



**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION  
DECISION SUMMARY  
ASSAY ONLY**

**I Background Information:**

**A 510(k) Number**

K233019

**B Applicant**

Hangzhou AllTest Biotech Co.,Ltd

**C Proprietary and Established Names**

AllTest Multi-Drug Rapid Test Cup; AllTest Multi-Drug Rapid Test Panel; AllTest Multi-Drug Rapid Test Cup Rx; AllTest Multi-Drug Rapid Test Panel Rx

**D Regulatory Information**

<b>Product Code(s)</b>	<b>Classification</b>	<b>Regulation Section</b>	<b>Panel</b>
DKZ	Class II	21 CFR 862.3100 - Amphetamine Test System	TX - Clinical Toxicology
NFT	Class II	21 CFR 862.3100 - Amphetamine test system	TX - Clinical Toxicology
DIS	Class II	21 CFR 862.3150 - Barbiturate test system	TX - Clinical Toxicology
JXM	Class II	21 CFR 862.3170 - Benzodiazepine test system	TX - Clinical Toxicology
DJG	Class II	21 CFR 862.3650 - Opiate test system	TX - Clinical Toxicology
PTH	Class II	21 CFR 862.3150 - Barbiturate test system	TX - Clinical Toxicology
DIO	Class II	21 CFR 862.3250 - Cocaine and cocaine metabolite test system	TX - Clinical Toxicology

DJR	Class II	21 CFR 862.3620 - Methadone test system	TX - Clinical Toxicology
NFV	Class II	21 CFR 862.3170 - Benzodiazepine test system	TX - Clinical Toxicology
LAF	Class II	21 CFR 862.3610 - Methamphetamine test system	TX - Clinical Toxicology
DJC	Class II	21 CFR 862.3610 - Methamphetamine test system	TX - Clinical Toxicology
LDJ	Class II	21 CFR 862.3870 - Cannabinoid test system	TX - Clinical Toxicology
NGL	Class II	21 CFR 862.3650 - Opiate test system	TX - Clinical Toxicology
NFY	Class II	21 CFR 862.3250 - Cocaine and cocaine metabolite test system	TX - Clinical Toxicology
DNK	Class II	21 CFR 862.3640 - Morphine test system	TX - Clinical Toxicology
LCM	Unclassified		
PTG	Class II	21 CFR 862.3620 - Methadone test system	TX - Clinical Toxicology
LFG	Class II	21 CFR 862.3910 - Tricyclic antidepressant drugs test system	TX - Clinical Toxicology
NGG	Class II	21 CFR 862.3610 - Methamphetamine test system	TX - Clinical Toxicology
NFW	Class II	21 CFR 862.3870 - Cannabinoid test system	TX - Clinical Toxicology
NGI	Class II	21 CFR 862.3640 - Morphine test system	TX - Clinical Toxicology
NGM	Unclassified		
QAW	Class II	21 CFR 862.3910 - Tricyclic antidepressant drugs test system	TX - Clinical Toxicology

## II Submission/Device Overview:

### A Purpose for Submission:

New Device

**B Measurand:**

Amphetamine, Secobarbital, Benzodiazepines, Buprenorphine, Cocaine, Marijuana, Methadone, Methamphetamine, Methylenedioxymethamphetamine, Morphine, Phencyclidine, Nortriptyline, Oxycodone, and 2-ethylidene-1,5-dimethyl-3,3- diphenylpyrrolidine.

**C Type of Test:**

Qualitative lateral flow immunoassay

**III Intended Use/Indications for Use:**

**A Intended Use(s):**

See Indications for Use below.

**B Indication(s) for Use:**

AllTest Multi-Drug Rapid Test Cup tests are competitive binding, lateral flow immunochromatographic assays for qualitative and simultaneous detection of Amphetamine, Buprenorphine, Secobarbital, Benzodiazepines, Cocaine, 2- ethylidene-1,5-dimethyl-3,3- diphenylpyrrolidine, Methamphetamine, Methylenedioxymethamphetamine, Morphine, Methadone, Oxycodone, Phencyclidine, Nortriptyline and Marijuana in human urine at the cutoff concentrations of:

<b>Drug (Identifier)</b>	<b>Cut-off level</b>
Amphetamine (AMP)	500 or 1000 ng/mL
Buprenorphine (BUP)	10 ng/mL
Secobarbital (BAR)	300 ng/mL
Benzodiazepines (BZO)	300 ng/mL
Cocaine (COC)	150 or 300 ng/mL
2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	300 ng/mL
Methamphetamine (MET)	500 or 1000 ng/mL
Methylenedioxymethamphetamine (MDMA)	500 ng/mL
Morphine (MOP/OPI)	300 or 2000 ng/mL
Methadone (MTD)	300 ng/mL
Oxycodone (OXY)	100 ng/mL
Phencyclidine (PCP)	25 ng/mL
Nortriptyline (TCA)	1000 ng/mL
Marijuana (THC)	50 ng/mL

AllTest Multi-Drug Rapid Test Cup offers any combinations of the above listed analytes. It is for in vitro diagnostic use only. It is intended for OTC use.

The tests may yield positive results for the prescription drugs Buprenorphine, Nortriptyline, Benzodiazepines, Secobarbital, and Oxycodone when taken at or above prescribed doses. It is not intended to distinguish between prescription use or abuse of these drugs. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly in evaluating a preliminary positive result.

The tests provide only preliminary results. To obtain a confirmed analytical result, a more specific alternate chemical method must be used. GC/MS or LC/MS is the recommended confirmatory method.

AllTest Multi-Drug Rapid Test Cup Rx tests are competitive binding, lateral flow immunochromatographic assays for qualitative and simultaneous detection of Amphetamine, Buprenorphine, Secobarbital, Benzodiazepines, Cocaine, 2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine, Methamphetamine, Methylenedioxymethamphetamine, Morphine, Methadone, Oxycodone, Phencyclidine, Nortriptyline and Marijuana in human urine at the cutoff concentrations of:

<b>Drug (Identifier)</b>	<b>Calibrator</b>	<b>Cut-off (ng/mL)</b>
Amphetamine (AMP)	d-Amphetamine	500 or 1000
Buprenorphine (BUP)	Buprenorphine	10
Secobarbital (BAR)	Secobarbital	300
Benzodiazepines (BZO)	Oxazepam	300
Cocaine (COC)	Benzoyllecgonine	150 or 300
2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine	300
Methamphetamine (MET)	d-Methamphetamine	500 or 1000
Methylenedioxymethamphetamine (MDMA)	d,l-Methylenedioxymethamphetamine	500
Morphine (MOP/OPI)	Morphine	300 or 2000
Methadone (MTD)	Methadone	300
Oxycodone (OXY)	Oxycodone	100
Phencyclidine (PCP)	Phencyclidine	25
Nortriptyline (TCA)	Nortriptyline	1000
Marijuana (THC)	11-nor- $\Delta^9$ -THC-9 COOH	50

AllTest Multi-Drug Rapid Test Cup Rx offers any combinations of the above listed analytes. It is for in vitro diagnostic use only. It is intended for prescription use.

The tests may yield positive results for the prescription drugs Buprenorphine, Nortriptyline, Benzodiazepines, Secobarbital, and Oxycodone when taken at or above prescribed doses. It is not intended to distinguish between prescription use or abuse of these drugs. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly in evaluating a preliminary positive result.

The tests provide only preliminary results. To obtain a confirmed analytical result, a more specific alternate chemical method must be used. GC/MS or LC/MS is the recommended confirmatory method.

AllTest Multi-Drug Rapid Test Panel tests are competitive binding, lateral flow immunochromatographic assays for qualitative and simultaneous detection of Amphetamine, Buprenorphine, Secobarbital, Benzodiazepines, Cocaine, 2- ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine, Methamphetamine, Methylenedioxymethamphetamine, Morphine, Methadone, Oxycodone, Phencyclidine, Nortriptyline and Marijuana in human urine at the cutoff concentrations of:

<b>Drug (Identifier)</b>	<b>Cut-off level</b>
Amphetamine (AMP)	500 or 1000 ng/mL
Buprenorphine (BUP)	10 ng/mL
Secobarbital (BAR)	300 ng/mL
Benzodiazepines (BZO)	300 ng/mL
Cocaine (COC)	150 or 300 ng/mL
2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	300 ng/mL
Methamphetamine (MET)	500 or 1000 ng/mL
Methylenedioxymethamphetamine (MDMA)	500 ng/mL
Morphine (MOP/OPI)	300 or 2000 ng/mL
Methadone (MTD)	300 ng/mL
Oxycodone (OXY)	100 ng/mL
Phencyclidine (PCP)	25 ng/mL
Nortriptyline (TCA)	1000 ng/mL
Marijuana (THC)	50 ng/mL

AllTest Multi-Drug Rapid Test Panel offers any combinations of the above listed analytes. It is for in vitro diagnostic use only. It is intended for OTC use.

The tests may yield positive results for the prescription drugs Buprenorphine, Nortriptyline, Benzodiazepines, Secobarbital, and Oxycodone when taken at or above prescribed doses. It is not intended to distinguish between prescription use or abuse of these drugs. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly in evaluating a preliminary positive result.

The tests provide only preliminary results. To obtain a confirmed analytical result, a more specific alternate chemical method must be used. GC/MS or LC/MS is the recommended confirmatory method.

AllTest Multi-Drug Rapid Test Panel Rx tests are competitive binding, lateral flow immunochromatographic assays for qualitative and simultaneous detection of Amphetamine, Buprenorphine, Secobarbital, Benzodiazepines, Cocaine, 2- ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine, Methamphetamine, Methylenedioxymethamphetamine, Morphine, Methadone, Oxycodone, Phencyclidine, Nortriptyline and Marijuana in human urine at the cutoff concentrations of:

<b>Drug (Identifier)</b>	<b>Calibrator</b>	<b>Cut-off (ng/mL)</b>
Amphetamine (AMP)	d-Amphetamine	500 or 1000

Buprenorphine (BUP)	Buprenorphine	10
Secobarbital (BAR)	Secobarbital	300
Benzodiazepines (BZO)	Oxazepam	300
Cocaine (COC)	Benzoyllecgonine	150 or 300
2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine	300
Methamphetamine (MET)	d-Methamphetamine	500 or 1000
Methylenedioxyamphetamine (MDMA)	d,1-Methylenedioxyamphetamine	500
Morphine (MOP/OPI)	Morphine	300 or 2000
Methadone (MTD)	Methadone	300
Oxycodone (OXY)	Oxycodone	100
Phencyclidine (PCP)	Phencyclidine	25
Nortriptyline (TCA)	Nortriptyline	1000
Marijuana (THC)	11-nor- $\Delta^9$ -THC-9 COOH	50

AllTest Multi-Drug Rapid Test Panel Rx offers any combinations of the above listed analytes. It is for in vitro diagnostic use only. It is intended for prescription use.

The tests may yield positive results for the prescription drugs Buprenorphine, Nortriptyline, Benzodiazepines, Secobarbital, and Oxycodone when taken at or above prescribed doses. It is not intended to distinguish between prescription use or abuse of these drugs. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly in evaluating a preliminary positive result.

The tests provide only preliminary results. To obtain a confirmed analytical result, a more specific alternate chemical method must be used. GC/MS or LC/MS is the recommended confirmatory method.

**C Special Conditions for Use Statement(s):**

Rx and OTC

**D Special Instrument Requirements:**

Not Applicable.

**IV Device/System Characteristics:**

**A Device Description:**

The AllTest Multi-Drug Rapid Test Cup and AllTest Multi-Drug Rapid Test Panel (Rx and OTC) are lateral flow-based, competitive immunochromatographic assays for the detection of Amphetamine, Secobarbital, Benzodiazepines, Buprenorphine, Cocaine, Marijuana, Methadone, Methamphetamine, Methylenedioxyamphetamine, Morphine, Phencyclidine, Nortriptyline, Oxycodone, and 2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine. The devices are intended for use as the first step in a two-step process to provide consumers with information concerning

the presence or absence of the above stated drugs in a urine sample. Information regarding confirmatory testing – the second step in the process is described in the instructions for use.

## **REAGENTS AND MATERIALS SUPPLIED**

(25 Tests per Kit)

- 25 Test Cups or Panels with desiccant in a sealed foil pouch
- 1 Instruction for use

## **MATERIALS REQUIRED BUT NOT PROVIDED**

- Timer
- Specimen Collection Container

### **B Principle of Operation:**

The AllTest Multi-Drug Rapid Test Cup/Panel tests are rapid tests for the qualitative detection of target drug or drug metabolites in urine samples. The tests are lateral flow chromatographic immunoassays. During testing, a urine specimen migrates upward by capillary action. If target drugs present in the urine specimen are below the cutoff concentration, it will not saturate the binding sites of its specific monoclonal mouse antibody coated on the particles. The antibody-coated particles will then be captured by immobilized drug-conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the target drug level exceeds its cutoff-concentration because it will saturate all the binding sites of the antibody coated on the particles. A band should form in the control region of the devices regardless of the presence of drug or metabolite in the sample to indicate that the tests have been performed properly.

### **V Substantial Equivalence Information:**

#### **A Predicate Device Name(s):**

Single and Multi-Drug Rapid Test Panel With Adulteration (Urine) Single and Multi-Drug Rapid Test Panel (Urine) Single and Multi-Drug Rapid Test Cup With Adulteration (Urine) Single and Multi-Drug Rapid Test Cup (Urine) Single Drug Rapid Test Dipstick (Urine) Single and Multi-Drug Home Rapid Test Panel (Urine) Single and Multi-Drug Home Rapid Test Cup (Urine) Single Drug Home Rapid Test Dipstick (Urine)

#### **B Predicate 510(k) Number(s):**

K182738

**C Comparison with Predicate(s):**

<b>Device &amp; Predicate Device(s):</b>	<u>K233019</u>	<u>K182738</u>
Device Trade Name	AllTest Multi-Drug Rapid Test Cup AllTest Multi-Drug Rapid Test Cup Rx AllTest Multi-Drug Rapid Test Panel AllTest Multi-Drug Rapid Test Panel Rx	Single and Multi-Drug Rapid Test Panel With Adulteration (Urine) Single and Multi-Drug Rapid Test Panel (Urine) Single and Multi-Drug Rapid Test Cup With Adulteration (Urine) Single and Multi-Drug Rapid Test Cup (Urine) Single Drug Rapid Test Dipstick (Urine) Single and Multi-Drug Home Rapid Test Panel (Urine) Single and Multi-Drug Home Rapid Test Cup (Urine) Single Drug Home Rapid Test Dipstick (Urine)
<b>General Device Characteristic Similarities</b>		
Intended Use/Indications For Use	For the qualitative determination of drugs and drug metabolites in human urine.	Same
Specimen Type	Human Urine	Same
Methodology	Competitive binding, lateral flow immunochromatographic assays based on the principle of antigen antibody immunochemistry.	Same
Configurations	Panel and Cup	Same
Intended Use	Rx and OTC	Same
<b>General Device Characteristic Differences</b>		

<p><b>Calibrator and Cutoff Values</b></p>	<p>Amphetamine (AMP): 500 or 1000 ng/ml Benzodiazepines (BZO):300 ng/ml Cocaine (COC): 150 or 300 ng/ml 11-Nor-<math>\Delta</math>9- Tetrahydrocannabinol-9- COOH (THC):50 ng/ml Methamphetamine (MET): 500 or 1000 ng/ml Morphine (MOP/OPI): 300 or 2000 ng/mL Oxycodone(OXY) : 100 ng/ml Secobarbital (BAR): 300 ng/ml Methadone (MTD): 300 ng/ml Buprenorphine (BUP): 10 ng/ml D,L- Methylenedioxymethamp hetamine (MDMA): 500 ng/ml Phencyclidine (PCP): 25 ng/ml Nortriptyline (TCA): 1000 ng/ml 2-ethylidene-1,5- dimethyl-3,3- diphenylpyrrolidine (EDDP): 300 ng/ml</p>	<p>Same except AMP 500 ng/ml only MET 500 ng/ml only COC 150 ng/ml only</p>
--	---	---

**VI Standards/Guidance Documents Referenced:**

None.

**VII Performance Characteristics (if/when applicable):**

**A Analytical Performance:**

1. Precision/Reproducibility:

Precision studies were carried out for samples with concentrations of -100% cutoff, -75% cutoff, -50% cutoff, -25% cutoff, cutoff, +25% cutoff, +50% cutoff, +75% cutoff and +100% cutoff using 3 lots of the test cup and 3 lots of the test panel. These samples were prepared by spiking drug in negative urine samples. Each drug concentration was confirmed by LC/MS. All sample aliquots were blindly labeled by the person who prepared the samples and didn't take part in the sample testing. For each concentration, tests were performed 5 replicates per day for 5 days per device in a randomized order. The results obtained for Amphetamine 1000, Cocaine 300, and Methamphetamine 1000 are summarized in the following tables. The rest of the data for Buprenorphine, Methylenedioxymethamphetamine, Secobarbital, Benzodiazepines, EDDP, Morphine, Methadone, Oxycodone, Phencyclidine, Nortriptyline and Marijuana were reported in K182738.

**AMP 1000 Cup**

Concen. by LC/MS (ng/mL)	-100% Cutoff	-75% Cutoff	-50% Cutoff	-25% Cutoff	Cutoff	Cutoff +25%	Cutoff +50%	Cutoff +75%	Cutoff +100%
Lot Number									
Lot 1	25-/0+	25-/0+	25-/0+	25-/0+	19+/6-	25+/0-	25+/0-	25+/0-	25+/0-
Lot 2	25-/0+	25-/0+	25-/0+	25-/0+	19+/6-	25+/0-	25+/0-	25+/0-	25+/0-
Lot 3	25-/0+	25-/0+	25-/0+	25-/0+	20+/5-	25+/0-	25+/0-	25+/0-	25+/0-

**AMP 1000 Panel**

Concen. by LC/MS (ng/mL)	-100% Cutoff	-75% Cutoff	-50% Cutoff	-25% Cutoff	Cutoff	Cutoff +25%	Cutoff +50%	Cutoff +75%	Cutoff +100%
Lot Number									
Lot 1	25-/0+	25-/0+	25-/0+	25-/0+	20+/5-	25+/0-	25+/0-	25+/0-	25+/0-
Lot 2	25-/0+	25-/0+	25-/0+	25-/0+	19+/6-	25+/0-	25+/0-	25+/0-	25+/0-
Lot 3	25-/0+	25-/0+	25-/0+	25-/0+	21+/4-	25+/0-	25+/0-	25+/0-	25+/0-

**COC 300 Cup**

Concen. by LC/MS (ng/mL)	-100% Cutoff	-75% Cutoff	-50% Cutoff	-25% Cutoff	Cutoff	Cutoff +25%	Cutoff +50%	Cutoff +75%	Cutoff +100%
Lot Number									

<b>Lot Number</b>									
<b>Lot 1</b>	25-/0+	25-/0+	25-/0+	25-/0+	20+/5-	25+/0-	25+/0-	25+/0-	25+/0-
<b>Lot 2</b>	25-/0+	25-/0+	25-/0+	25-/0+	19+/6-	25+/0-	25+/0-	25+/0-	25+/0-
<b>Lot 3</b>	25-/0+	25-/0+	25-/0+	25-/0+	21+/4-	25+/0-	25+/0-	25+/0-	25+/0-

**COC 300 Panel**

<b>Concen. by LC/MS (ng/mL)</b>	<b>-100% Cutoff</b>	<b>-75% Cutoff</b>	<b>-50% Cutoff</b>	<b>-25% Cutoff</b>	<b>Cutoff</b>	<b>Cutoff +25%</b>	<b>Cutoff +50%</b>	<b>Cutoff +75%</b>	<b>Cutoff +100%</b>
<b>Lot Number</b>									
<b>Lot 1</b>	25-/0+	25-/0+	25-/0+	25-/0+	20+/5-	25+/0-	25+/0-	25+/0-	25+/0-
<b>Lot 2</b>	25-/0+	25-/0+	25-/0+	25-/0+	20+/5-	25+/0-	25+/0-	25+/0-	25+/0-
<b>Lot 3</b>	25-/0+	25-/0+	25-/0+	25-/0+	20+/5-	25+/0-	25+/0-	25+/0-	25+/0-

**MET 1000 Cup**

<b>Concen. by LC/MS (ng/mL)</b>	<b>-100% Cutoff</b>	<b>-75% Cutoff</b>	<b>-50% Cutoff</b>	<b>-25% Cutoff</b>	<b>Cutoff</b>	<b>Cutoff +25%</b>	<b>Cutoff +50%</b>	<b>Cutoff +75%</b>	<b>Cutoff +100%</b>
<b>Lot Number</b>									
<b>Lot 1</b>	25-/0+	25-/0+	25-/0+	25-/0+	19+/6-	25+/0-	25+/0-	25+/0-	25+/0-
<b>Lot 2</b>	25-/0+	25-/0+	25-/0+	25-/0+	19+/6-	25+/0-	25+/0-	25+/0-	25+/0-
<b>Lot 3</b>	25-/0+	25-/0+	25-/0+	25-/0+	20+/5-	25+/0-	25+/0-	25+/0-	25+/0-

**MET 1000 Panel**

<b>Concen. by LC/MS (ng/mL)</b>	<b>-100% Cutoff</b>	<b>-75% Cutoff</b>	<b>-50% Cutoff</b>	<b>-25% Cutoff</b>	<b>Cutoff</b>	<b>Cutoff +25%</b>	<b>Cutoff +50%</b>	<b>Cutoff +75%</b>	<b>Cutoff +100%</b>
<b>Lot Number</b>									

<b>Lot 1</b>	25-/0+	25-/0+	25-/0+	25-/0+	20+/5-	25+/0-	25+/0-	25+/0-	25+/0-
<b>Lot 2</b>	25-/0+	25-/0+	25-/0+	25-/0+	21+/4-	25+/0-	25+/0-	25+/0-	25+/0-
<b>Lot 3</b>	25-/0+	25-/0+	25-/0+	25-/0+	21+/4-	25+/0-	25+/0-	25+/0-	25+/0-

2. Linearity:

Not Applicable.

3. Analytical Specificity/Interference:

**Interfering Substances**

Potential interfering substances found in human urine of physiological or pathological conditions were added to drug-free urine and target drugs urine with concentrations at 25% below and 25% above all analyte cut-off levels in the test. These urine samples were tested using three lots of each device format by twelve operators (six for the cup device and six for the panel device). Compounds that showed no interference at a concentration of 100µg/mL are summarized in the following tables. No differences were observed for the different device formats.

Acetylsalicylic Acid	5, 5-Diphenylhydantoin	19-Norethindrone
Albumin (100mg/dL)	Erythromycin	Noscipine
Amoxicillin	Estradiol	Octopamine
Ampicillin	Estrone	Papaverine
Aspartame	Ethanol (1%)	Penicillin-G
Aspirin	Fenofibrate	Perphenazine
Atropine	Fentanyl	Phenelzine
Baclofen	Fotemustine	Phenylethylamine
Benzocaine	Furosemide	Promazine
Benzoic Acid	Gemfibrozil	Promethazine
Bilirubin	Gentisic acid	Pyridoxine
Carisoprodol	Glucose	Pyrilamine
Chloramphenicol	Guaiacol glyceryl ether	Pyrogallol
Chlordiazepoxide	Hemoglobin	Quinine
(+)-Chlorpheniramine	Hydralazine	Quinolinic Acid
Chlorpromazine	Hydrocortisone	R-(-)-Apomorphine
Cholesterol	3-Hydroxytyramine	Ranitidine
Clonidine	(+/-)-Isoproterenol	Salicylic Acid
Cortisone	Ketamine	Sulindac
(-)-Cotinine	L-Ascorbic Acid	Tetracycline
Creatine Hydrate	Meprobamate	Tetrahydrozoline
Creatinine	Methylphenidate	Thiamine
Cyclodextrin	Nalidixic Acid	Thioridazine
d,l-Propranolol	Naltrexone	Tramadol
Deoxycorticosterone	(+)-Naproxen	Trifluoperazine
Dextromethorphan	Niacinamide	Tryptamine
Diclofenac	Nicotinic Acid	Uric Acid

**Cross Reactivity**

To test specificity, drug metabolites and other structurally related compounds that are likely to cross-react in urine samples were spiked into negative urine and were tested using three lots of each device. The lowest concentrations that caused a positive result for each compound are listed below for Amphetamine 1000, Cocaine 300, and Methamphetamine 1000. No differences were observed for different device format.

The data for Buprenorphine, Methylenedioxyamphetamine, Secobarbital, Benzodiazepines, Morphine, Methadone, Oxycodone, Phencyclidine, EDDP, Nortriptyline and Marijuana were reported in K182738.

<b>AMP 1000 (Cutoff=1000 ng/mL)</b>	<b>Result Positive at (ng/ml)</b>	<b>%Cross-Reactivity</b>
d-Amphetamine	1000	100%
Methylenedioxyethylamphetamine (MDEA)	>100000	<1%
d,l-Methamphetamine	>100000	<1%
Phenylephrine	>100000	<1%
d-Methamphetamine	>100000	<1%
l-Methamphetamine	>100000	<1%
d,l - Methylenedioxy methamphetamine	>100000	<1%
l-Amphetamine	1000	100%
Ephedrine	>100000	<1%
Pseudoephedrine	>100000	<1%
d, l-Amphetamine	1000	100%
d,l-3,4-Methylenedioxyamphetamine (MDA)	50	2000%
Phentermine	1000	100%

<b>COC 300 (Cutoff=300 ng/mL)</b>	<b>Result Positive at (ng/ml)</b>	<b>%Cross-Reactivity</b>
Benzoyllecgonine	300	100%
Cocaine	250	120%
Cocaethylene	500	60%
Ecgonine	>100000	<0.3%
Norcocaine	>100000	<0.3%

<b>MET 1000 (Cutoff=1000 ng/mL)</b>	<b>Result Positive at (ng/ml)</b>	<b>%Cross-Reactivity</b>
d-Methamphetamine	1000	100%
l -Methamphetamine	25000	4%
d,l-Amphetamine	500	200%

Phentermine	>100000	<1%
d,l-Methamphetamine	500	200%
d-Amphetamine	>100000	<1%
l-Amphetamine	>100000	<1%
Ephedrine	>100000	<1%
Phenylephrine	>100000	<1%
Pseudoephedrine	>100000	<1%
3,4- Methylendioxy methamphetamine (MDMA)	2500	40%
d,l- Methylendioxy ethylamphetamine (MDEA)	12500	8.0%
d,l-3,4-Methylendioxyamphetamine (MDA)	>100000	<1%

### **Effect of Urine Specific Gravity and Urine pH**

Effect of urine density and pH data were reported only for Amphetamine 1000, Cocaine 300 and Methamphetamine 1000. The rest of the data were reported in the cleared AllTest submission of K182738.

The effect of pH was evaluated by preparing samples with test drug analytes at +/-25% of the cutoff for all devices and adjusting the pH to 4.00 to 9.00, in increments of 1 pH unit. There were no interferences observed.

The effect of specific gravity was evaluated by preparing samples with test drug analytes at +/-25% of the cutoff for all devices and adjusting the specific gravity to 1.000, 1.005, 1.008, 1.011, 1.013, 1.016, 1.022, 1.025, 1.027, 1.029, 1.033, 1.035. There were no interferences observed.

4. Assay Reportable Range:

Not Applicable.

5. Traceability, Stability, Expected Values (Controls, Calibrators, or Methods):

All drug calibrators of the device are traceable to available commercial reference materials.

6. Detection Limit:

See precision data in Section VII.A.1., above, for assay performance around the claimed cutoff concentrations.

7. Assay Cutoff:

Characterization of how the device performs analytically around the claimed cutoff concentration appears in the precision section VII.A.1., above.

**B Comparison Studies:**

1. Method Comparison with Predicate Device:

Method comparison studies were performed at the manufacturer’s site using three operators for each device format (cup and panel). Operators tested a minimum of 40 unaltered positive and 40 unaltered negative clinical samples for each drug. The samples were blind labeled and compared to LC/MS results.

The results are presented in the tables below for Amphetamine 1000, Cocaine 300, and Methamphetamine 1000. The data for Buprenorphine, Methylenedioxymethamphetamine, Secobarbital, Benzodiazepines, Morphine, Methadone, Oxycodone, Phencyclidine, EDDP, Nortriptyline and Marijuana were reported in K182738.

**AMP Cup**

AllTest Multi-Drug Test Cup		Negative	Low Negative by LC/MS (less than -50%)	Near Cutoff Negative by LC/MS (Between -50% and cutoff)	Near Cutoff Positive by LC/MS (Between the cutoff and +50%)	High Positive by LC/MS (greater than +50%)
Viewer A	Positive	0	0	1	13	26
	Negative	12	12	15	1	0
Viewer B	Positive	0	0	0	13	26
	Negative	12	12	16	1	0
Viewer C	Positive	0	0	1	14	26
	Negative	12	12	15	0	0

**Discordant Results**

Viewer	Sample Number	LC/MS Result	Viewer Results
Viewer A	SN165	896.471	+
Viewer A	SN175	1070.815	-
Viewer B	SN044	1149.522	-
Viewer C	SN128	838.956	+

### AMP Panel

AllTest Multi-Drug Test Panel		Negative	Low Negative by LC/MS (less than -50%)	Near Cutoff Negative by LC/MS (Between -50% and cutoff)	Near Cutoff Positive by LC/MS (Between the cutoff and +50%)	High Positive by LC/MS (greater than +50%)
Viewer A	Positive	0	0	0	13	26
	Negative	12	12	16	1	0
Viewer B	Positive	0	0	1	14	26
	Negative	12	12	15	0	0
Viewer C	Positive	0	0	2	13	26
	Negative	12	12	14	1	0

### Discordant Results

Viewer	Sample Number	LC/MS Result	Viewer Results
Viewer A	SN044	1149.522	-
Viewer B	SN062	833.996	+
Viewer C	SN036	922.995	+
Viewer C	SN120	987.795	+
Viewer C	SN150	1181.731	-

### COC Cup

AllTest Multi-Drug Test Cup		Negative	Low Negative by LC/MS (less than -50%)	Near Cutoff Negative by LC/MS (Between -50% and cutoff)	Near Cutoff Positive by LC/MS (Between the cutoff and +50%)	High Positive by LC/MS (greater than +50%)
Viewer A	Positive	0	0	0	9	30
	Negative	12	10	18	1	0
Viewer B	Positive	0	0	1	10	30
	Negative	12	10	17	0	0
Viewer C	Positive	0	0	1	10	30
	Negative	12	10	17	0	0

Discordant Results

Viewer	Sample Number	LC/MS Result	Viewer Results
Viewer A	SN174	324.6	-
Viewer B	SN070	297.0	+
Viewer C	SN079	297.7	+

**COC Panel**

AllTest Multi-Drug Test Panel		Negative	Low Negative by LC/MS (less than -50%)	Near Cutoff Negative by LC/MS (Between -50% and cutoff)	Near Cutoff Positive by LC/MS (Between the cutoff and +50%)	High Positive by LC/MS (greater than +50%)
Viewer A	Positive	0	0	1	8	30
	Negative	12	10	17	2	0
Viewer B	Positive	0	0	1	10	30
	Negative	12	10	17	0	0
Viewer C	Positive	0	0	0	8	30
	Negative	12	10	18	2	0

Discordant Results

Viewer	Sample Number	LC/MS Result	Viewer Results
Viewer A	SN079	297.79	+
Viewer A	SN142	352.11	-
Viewer A	SN174	324.69	-
Viewer B	SN012	292.74	+
Viewer C	SN084	312.95	-
Viewer C	SN142	352.11	-

**MET Cup**

AllTest Multi-Drug Test Cup		Negative	Low Negative by LC/MS (less than -50%)	Near Cutoff Negative by LC/MS (Between -50% and cutoff)	Near Cutoff Positive by LC/MS (Between the cutoff and +50%)	High Positive by LC/MS (greater than +50%)
Viewer A	Positive	0	0	1	10	29
	Negative	12	11	16	1	0
Viewer B	Positive	0	0	0	10	29

AllTest Multi-Drug Test Cup		Negative	Low Negative by LC/MS (less than -50%)	Near Cutoff Negative by LC/MS (Between -50% and cutoff)	Near Cutoff Positive by LC/MS (Between the cutoff and +50%)	High Positive by LC/MS (greater than +50%)
	Negative	12	11	17	1	0
Viewer C	Positive	0	0	2	10	29
	Negative	12	11	15	1	0

Discordant Results

Viewer	Sample Number	LC/MS Result	Viewer Results
Viewer A	SN108	974.577	+
Viewer A	SN043	1194.980	-
Viewer B	SN043	1194.980	-
Viewer C	SN054	919.834	+
Viewer C	SN188	841.261	+
Viewer C	SN197	1167.315	-

MET Panel

AllTest Multi-Drug Test Panel		Negative	Low Negative by LC/MS (less than -50%)	Near Cutoff Negative by LC/MS (Between -50% and cutoff)	Near Cutoff Positive by LC/MS (Between the cutoff and +50%)	High Positive by LC/MS (greater than +50%)
Viewer A	Positive	0	0	1	9	29
	Negative	12	11	16	2	0
Viewer B	Positive	0	0	1	10	29
	Negative	12	11	16	1	0
Viewer C	Positive	0	0	0	10	29
	Negative	12	11	17	1	0

Discordant Results

Viewer	Sample Number	LC/MS Result	Viewer Results
Viewer A	SN159	882.244	+
Viewer A	SN043	1194.980	-
Viewer A	SN186	1138.254	-
Viewer B	SN159	882.244	+

<b>Viewer B</b>	SN186	1138.254	-
<b>Viewer C</b>	SN197	1167.315	-

### Lay User Study

A lay user study was performed at three intended user sites by 560 lay users. Each lay user tested only one sample (140 tested the low cutoff test cup, 140 tested the low cutoff test panel, 140 tested the high cutoff test cup, and 140 tested the high cutoff test panel). The lay users had diverse educational and professional backgrounds and ranged in age from 20 to > 50 years. Urine samples were prepared at the following concentrations; negative, +/-75%, +/-50%, +/-25% of the cutoff by spiking drugs into drug free-pooled urine specimens. The concentrations of the samples were confirmed by LC/MS. Each sample was aliquoted into individual containers and blind-labeled. Each participant was provided with the package insert, 1 blind labeled sample and a device. Each device format was tested. Results are shown below.

### **Results for Low Cutoff Cup**

<b>Drug</b>	<b>Cutoff (ng/mL)</b>	<b>Results</b>	<b>Concentration</b>						
			<b>-100% cutoff</b>	<b>-75% cutoff</b>	<b>-50% cutoff</b>	<b>-25% cutoff</b>	<b>+25% cutoff</b>	<b>+50% cutoff</b>	<b>+75% cutoff</b>
AMP	500	Negative	20	20	20	19	2	0	0
		Positive	0	0	0	1	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	90%	100%	100%
BAR	300	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
BZO	300	Negative	20	20	20	18	1	0	0
		Positive	0	0	0	2	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	95%	100%	100%
BUP	10	Negative	20	20	20	18	1	0	0
		Positive	0	0	0	2	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	95%	100%	100%
COC	150	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
		Negative	20	20	20	18	2	0	0

Drug	Cutoff (ng/mL)	Results	Concentration						
			-100% cutoff	-75% cutoff	-50% cutoff	-25% cutoff	+25% cutoff	+50% cutoff	+75% cutoff
EDDP	300	Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
		Negative	20	20	20	18	1	0	0
MDMA	500	Positive	0	0	0	2	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	95%	100%	100%
		Negative	20	20	20	18	2	0	0
MET	500	Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
		Negative	20	20	20	18	2	0	0
MOP	300	Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
		Negative	20	20	20	18	2	0	0
MTD	300	Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
		Negative	20	20	20	18	2	0	0
OXY	100	Positive	0	0	0	1	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	90%	100%	100%
		Negative	20	20	20	19	2	0	0
PCP	25	Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
		Negative	20	20	20	18	2	0	0
TCA	1000	Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
		Negative	20	20	20	18	2	0	0
THC	50	Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
		Negative	20	20	20	18	2	0	0

### Results for High Cutoff Cup

Drug	Cutoff (ng/mL)	Results	Concentration						
			-100% cutoff	-75% cutoff	-50% cutoff	-25% cutoff	+25% cutoff	+50% cutoff	+75% cutoff
AMP	1000	Negative	20	20	20	18	1	0	0
		Positive	0	0	0	2	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	95%	100%	100%
BAR	300	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
BZO	300	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
BUP	10	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
COC	300	Negative	20	20	20	18	1	0	0
		Positive	0	0	0	2	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	95%	100%	100%
EDDP	300	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
MDMA	500	Negative	20	20	20	18	1	0	0
		Positive	0	0	0	2	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	95%	100%	100%
MET	1000	Negative	20	20	20	19	1	0	0
		Positive	0	0	0	1	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	95%	100%	100%

Drug	Cutoff (ng/mL)	Results	Concentration						
			-100% cutoff	-75% cutoff	-50% cutoff	-25% cutoff	+25% cutoff	+50% cutoff	+75% cutoff
OPI	2000	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
MTD	300	Negative	20	20	20	19	1	0	0
		Positive	0	0	0	1	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	95%	100%	100%
OXY	100	Negative	20	20	20	18	1	0	0
		Positive	0	0	0	2	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	95%	100%	100%
PCP	25	Negative	20	20	20	19	2	0	0
		Positive	0	0	0	1	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	90%	100%	100%
TCA	1000	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
THC	50	Negative	20	20	20	19	2	0	0
		Positive	0	0	0	1	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	90%	100%	100%

#### Results for Low Cutoff Panel

Drug	Cutoff (ng/mL)	Results	Concentration						
			-100% cutoff	-75% cutoff	-50% cutoff	-25% cutoff	+25% cutoff	+50% cutoff	+75% cutoff
AMP	500	Negative	20	20	20	19	2	0	0
		Positive	0	0	0	1	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	90%	100%	100%
BAR	300	Negative	20	20	20	18	1	0	0
		Positive	0	0	0	2	19	20	20
		Total	20	20	20	20	20	20	20

Drug	Cutoff (ng/mL)	Results	Concentration						
			-100% cutoff	-75% cutoff	-50% cutoff	-25% cutoff	+25% cutoff	+50% cutoff	+75% cutoff
		Agreement (%)	100%	100%	100%	90%	95%	100%	100%
BZO	300	Negative	20	20	20	18	1	0	0
		Positive	0	0	0	2	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	95%	100%	100%
BUP	10	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
COC	150	Negative	20	20	20	18	1	0	0
		Positive	0	0	0	2	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	95%	100%	100%
EDDP	300	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
MDMA	500	Negative	20	20	20	19	1	0	0
		Positive	0	0	0	1	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	95%	100%	100%
MET	500	Negative	20	20	20	18	1	0	0
		Positive	0	0	0	2	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	95%	100%	100%
MOP	300	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
MTD	300	Negative	20	20	20	18	1	0	0
		Positive	0	0	0	2	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	95%	100%	100%
		Negative	20	20	20	19	2	0	0

Drug	Cutoff (ng/mL)	Results	Concentration						
			-100% cutoff	-75% cutoff	-50% cutoff	-25% cutoff	+25% cutoff	+50% cutoff	+75% cutoff
OXY	100	Positive	0	0	0	1	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	90%	100%	100%
		Negative	20	20	20	19	2	0	0
PCP	25	Positive	0	0	0	1	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	90%	100%	100%
		Negative	20	20	20	19	2	0	0
TCA	1000	Positive	0	0	0	1	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	90%	100%	100%
		Negative	20	20	20	18	2	0	0
THC	50	Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
		Negative	20	20	20	18	2	0	0

#### Results for High Cutoff Panel

Drug	Cutoff (ng/mL)	Results	Concentration						
			-100% cutoff	-75% cutoff	-50% cutoff	-25% cutoff	+25% cutoff	+50% cutoff	+75% cutoff
AMP	1000	Positive	0	0	0	2	20	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	100%	100%	100%
		Negative	20	20	20	18	0	0	0
BAR	300	Positive	0	0	0	1	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	95%	100%	100%
		Negative	20	20	20	19	1	0	0
BZO	300	Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
		Negative	20	20	20	18	2	0	0
BUP	10	Positive	0	0	0	1	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	90%	100%	100%
		Negative	20	20	20	19	2	0	0

Drug	Cutoff (ng/mL)	Results	Concentration						
			-100% cutoff	-75% cutoff	-50% cutoff	-25% cutoff	+25% cutoff	+50% cutoff	+75% cutoff
COC	300	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
EDDP	300	Negative	20	20	20	19	2	0	0
		Positive	0	0	0	1	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	90%	100%	100%
MDMA	500	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
MET	1000	Negative	20	20	20	19	1	0	0
		Positive	0	0	0	1	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	95%	100%	100%
OPI	2000	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
MTD	300	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
OXY	100	Negative	20	20	20	19	1	0	0
		Positive	0	0	0	1	19	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	95%	100%	100%
PCP	25	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
TCA	1000	Negative	20	20	20	18	2	0	0
		Positive	0	0	0	2	18	20	20
		Total	20	20	20	20	20	20	20

Drug	Cutoff (ng/mL)	Results	Concentration						
			-100% cutoff	-75% cutoff	-50% cutoff	-25% cutoff	+25% cutoff	+50% cutoff	+75% cutoff
		Agreement (%)	100%	100%	100%	90%	90%	100%	100%
THC	50	Negative	20	20	20	19	2	0	0
		Positive	0	0	0	1	18	20	20
		Total	20	20	20	20	20	20	20
		Agreement (%)	100%	100%	100%	95%	90%	100%	100%

Lay users were also given surveys on the ease of understanding the package insert instructions. All lay users indicated that the device instructions can be easily followed. A Flesch-Kincaid reading analysis was performed on each package insert and the scores revealed a reading Grade Level of 7.

2. Matrix Comparison:

Not Applicable.

**C Clinical Studies:**

1. Clinical Sensitivity:

Not Applicable.

2. Clinical Specificity:

Not Applicable.

3. Other Clinical Supportive Data (When 1. and 2. Are Not Applicable):

Not Applicable.

**D Clinical Cutoff:**

Not Applicable.

**E Expected Values/Reference Range:**

Not Applicable.

**VIII Proposed Labeling:**

The labeling supports the finding of substantial equivalence for this device.

**IX Conclusion:**

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.