

K030534

MAY 20 2003

510(k) Summary

Sponsor: Skeletal Kinetics, LLC
10201 Bubb Road, Cupertino, CA 95014
Contact Person: Duran Yetkinler, M.D., Ph.D.
Phone Number: 408 366 5002
Fax Number: 408 366 1077
Prepared: February 18, 2003

Trade Name: Callos™
Common Name: Bone Graft Substitute
Classification: Unclassified
Product Code: 87 MVQ

Predicate Device: Callos Bone Void Filler is substantially equivalent to Norian SRS Bone Void Filler (K011897).

Device Description: Callos Bone Void Filler is an injectable, moldable and biocompatible bone void filler. Callos Bone Void Filler resorbs and is replaced with bone during the healing process. The 3 cc, 5 cc, and 10 cc Callos Bone Void Filler kits are provided sterile and are for single use only.

Intended Use/Indications for Use: Callos Bone Void Filler is indicated to fill bony voids or gaps of the skeletal system (i.e. extremities, spine, and pelvis). These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. Callos Bone Void Filler is indicated only for bony voids or gaps that are not intrinsic to the stability of the bony structure. The product provides a bone void filler that resorbs and is replaced by bone during the healing process.

Technological Characteristics: Similar to the predicate device, Callos Bone Void Filler is an injectable, moldable, biocompatible, resorbable calcium phosphate based material intended for identical indications.

Performance Data: Non-clinical testing included material properties such as density, porosity, dimensional stability, injectability, setting time, working time, pH, and setting temperature. Biocompatibility testing demonstrated that the material is non-cytotoxic, non-systemic toxic, non-mutagenic, non-irritative, non-pyrogenic, and non-sensitizing. Comparative testing with the predicate device showed equivalence in terms of solubility and dissolution rate, X-Ray Diffraction (XRD), Fourier Transform Infrared (FTIR) spectroscopy and elemental analysis. Animal testing demonstrated substantial equivalence to the predicate device following *in vivo* implantation. Histological, chemical, crystallographical, and mechanical analyses showed substantial equivalence.

Basis for Substantial Equivalence: The Callos Bone Void Filler has the same intended use, identical indications, and very similar technological characteristics as the predicate device. Any minor technological differences between Callos Bone Void Filler and its predicate device do not raise any new issues of safety or effectiveness.

Functional, biocompatibility, and animal testing results show that the Callos Bone Void Filler is as safe and effective as the predicate device. Thus, the Callos Bone Void Filler is substantially equivalent.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

MAY 20 2003

Skeletal Kinetics, LLC
c/o Mr. Howard Holstein
Regulatory Counsel
Hogan & Hartson, LLC
555 13th Street, N.W.
Washington, DC 20004

Re: K030554

Trade/Device Name: Callos[™] Bone Void Filler
Regulatory Class: Unclassified
Product Code: MQV
Dated: February 20, 2003
Received: February 21, 2003

Dear Mr. Holstein:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Mr. Howard Holstein

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read "Celia M. Witten". The signature is written in a cursive style with a large initial "C".

Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Indications for Use

Page 1 of 1

510(k) Number (if known): K030554

Device Name: Callos Bone Void Filler


Indications for use:

Callos Bone Void Filler is indicated to fill bony voids or gaps of the skeletal system (i.e. extremities, spine, pelvis). These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. Callos Bone Void Filler is indicated only for bony voids or gaps that are not intrinsic to the stability of the bony structure. The product provides a bone void filler that resorbs and is replaced by bone during the healing process.

(PLEASE DO NOT WRITE BELOW THIS LINE—CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use OR Over-the-Counter Use
(per 21 CFR 801.109)


(Division Sign-Off)
Division of General, Restorative and
Neurological Devices

510(k) Number K030554

Add to file



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JUN 24 2003

Duran N. Yetkinler, M.D., Ph.D.
VP of Regulatory and Product Development
Skeletal Kinetics, LLC
10201 Bubb Road
Cupertino, CA 95014

Re: K030554
Trade/Device Name: Callos™ Bone Void Filler
Regulatory Class: Unclassified
Product Code: MQV
Dated: February 20, 2003
Received: February 21, 2003

Dear Dr. Yetkinler:

This letter corrects our substantially equivalent letter of May 20, 2003, regarding the contact and address.

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent [(for the indications for use stated in the enclosure)] to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other


Page 2 – Dr. Duran N. Yetkinler

Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to continue marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at their toll free number (800) 638-2041 or at (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "for Mark N. Miller" with a flourish at the end.

Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Sloan, Nadine Y.

From: duran yetkinler [duran@skeletalkinetics.com]
Sent: Monday, June 09, 2003 4:31 PM
To: 'Sloan, Nadine Y.'
Subject: RE: Primary contact person and address change request

Nadine,

Hope all is well. I was following up on the request that I made regarding changing the name and address of the contact person. Any improvements? It turns out that we really want to make this change, so please advise me if I need to write any official letter. Thank you.

Best Regards,

Duran

A Duran N Yetkinler MD PhD
VP of Regulatory and Product Development
Skeletal Kinetics, LLC
10201 Bubb Road
Cupertino, CA 95014
Phone 408 366 5002
Fax 408 366 1077
Cell 408 757 6603

-----Original Message-----

From: Sloan, Nadine Y. [mailto:NYR@CDRH.FDA.GOV]
Sent: Friday, May 23, 2003 10:50 AM
To: 'duran yetkinler'
Subject: RE: Primary contact person and address change request

Duran,

Thank you for the information. I will look into this further and see about issuing a correction letter. I will get back to you sometime next week.

Sincerely,

Nadine

-----Original Message-----

From: duran yetkinler [mailto:duran@skeletalkinetics.com]
Sent: Friday, May 23, 2003 1:59 PM
To: NYR@CDRH.FDA.GOV
Subject: Primary contact person and address change request

Nadine,

I have received the substantial equivalence letter, K030554 dated May 20, 2003, for our product Callos. I noticed on the letter that the letter was addressed to Howard Holstein and his address. Howard Holstein and his group act as our regulatory consultant to us, but we would like to keep the primary contact person and address for Skeletal Kinetics as myself and Skeletal Kinetics address. The full name and address for Skeletal Kinetics are given below. I am requesting to make these changes on our files, so that any future and present confusions will be avoided. Thank you for your

6/9/2003

cooperation.

Sincerely,

Duran

Duran N Yetkinler MD PhD
VP of Regulatory and Product Development
Skeletal Kinetics, LLC
10201 Bubb Road
Cupertino, CA 95014
Phone 408 366 5002
Fax 408 366 1077
Cell 408 757 6603

6/9/2003

Dmc
JUL 14 2018



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

MAY 20 2003

Skeletal Kinetics, LLC
c/o Mr. Howard Holstein
Regulatory Counsel
Hogan & Hartson, LLC
555 13th Street, N.W.
Washington, DC 20004

Re: K030554

Trade/Device Name: Callos[™] Bone Void Filler
Regulatory Class: Unclassified
Product Code: MQV
Dated: February 20, 2003
Received: February 21, 2003

Dear Mr. Holstein:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Mr. Howard Holstein

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,


for Celia M. Witten, Ph.D., M.D.

Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Indications for Use

Page 1 of 1

510(k) Number (if known): K030554

Device Name: Callos Bone Void Filler


Indications for use:

Callos Bone Void Filler is indicated to fill bony voids or gaps of the skeletal system (i.e. extremities, spine, pelvis). These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. Callos Bone Void Filler is indicated only for bony voids or gaps that are not intrinsic to the stability of the bony structure. The product provides a bone void filler that resorbs and is replaced by bone during the healing process.

(PLEASE DO NOT WRITE BELOW THIS LINE—CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use OR Over-the-Counter Use
(per 21 CFR 801.109)


(Division Sign-Off)
Division of General, Restorative and
Neurological Devices

510(k) Number K030554

3

DEPARTMENT OF HEALTH AND HUMAN SERVICES

Public Health Service

Food and Drug Administration
Center for Devices and
Radiological Health
Office of Device Evaluation
Document Mail Center (HFZ-401)
9200 Corporate Blvd.
Rockville, Maryland 20850

February 21, 2003

SKELETAL KINETICS, LLC
C/O HOGAN & HARTSON, LLC
555 13TH STREET N.W.
WASHINGTON, DC 20004
ATTN: HOWARD HOLSTEIN

510(k) Number: K030554
Received: 21-FEB-2003
Product: CALLOS

The Center for Devices and Radiological Health (CDRH), Office of Device Evaluation (ODE), has received the Premarket Notification you submitted in accordance with Section 510(k) of the Federal Food, Drug, and Cosmetic Act (Act) for the above referenced product. We have assigned your submission a unique 510(k) number that is cited above. Please refer prominently to this 510(k) number in any future correspondence that relates to this submission. We will notify you when the processing of your premarket notification has been completed or if any additional information is required. YOU MAY NOT PLACE THIS DEVICE INTO COMMERCIAL DISTRIBUTION UNTIL YOU RECEIVE A LETTER FROM FDA ALLOWING YOU TO DO SO.

As a reminder, we would like to mention that FDA requires all 510(k) submitters to provide an indications for use statement on a separate page. If you have not included this indications for use statement in addition to your 510(k) summary (807.92), or a 510(k) statement (807.93), and your Truthful and Accurate statement, please do so as soon as possible. If the above mentioned requirements have been submitted, please do not submit them again. There may be other regulations or requirements affecting your device such as Postmarket Surveillance (Section 522(a)(1) of the Act) and the Device Tracking regulation (21 CFR Part 821). Please contact the Division of Small Manufacturers, International and Consumer Assistance (DSMICA) at the telephone or web site below for more information.

Please remember that all correspondence concerning your submission MUST be sent to the Document Mail Center (DMC)(HFZ-401) at the above letterhead address. Correspondence sent to any address other than the DMC will not be considered as part of your official premarket notification submission. Also, please note the new Blue Book Memorandum regarding Fax and E-mail Policy entitled, "Fax and E-Mail Communication with Industry about Premarket Files Under Review. Please refer to this guidance for information on current fax and e-mail practices at www.fda.gov/cdrh/ode/a02-01.html.

You should be familiar with the manual entitled, "Premarket Notification 510(k) Regulatory Requirements for Medical Devices" available from DSMICA. If you have other procedural or policy questions, or want information on how to check on the status of your submission (after 90 days from the receipt date), please contact DSMICA at (301) 443-6597 or its toll-free number (800) 638-2041, or at their Internet address <http://www.fda.gov/cdrh/dsmamain.html> or me at (301) 594-1190.

Sincerely yours,

Marjorie Shulman
Consumer Safety Officer
Premarket Notification Staff
Office of Device Evaluation
Center for Devices and Radiological Health

21

K 030554

February 20, 2003

Food and Drug Administration
Center for Devices and Radiological Health
Document Mail Center, HFZ-401
9200 Corporate Drive
Rockville, MD 20850

FDA/CDRH/REG-13
2003 FEB 21 A 11:59

Re: 510(k) Notification—Bone Void Filler

Dear Madam/Sir:

Enclosed is a submission pursuant to Section 510(k) of the Federal Food, Drug, and Cosmetic Act and the regulations contained in 21 CFR 807.

The purpose of this submission is to notify the FDA of our intent to market a bone void filler. Skeletal Kinetics' bone void filler, Callos™, is indicated to fill bony voids or gaps of the skeletal system (i.e. extremities, spine, pelvis). These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. Callos Bone Void Filler is indicated only for bony voids or gaps that are not intrinsic to the stability of the bony structure. The product provides a bone void filler that resorbs and is replaced by bone during the healing process. Callos will be marketed as a kit comprised of a powder (calcium phosphate source), liquid (dilute sodium silicate), mixing bowl, pestle and spatula. The user empties the contents of both vials into the mixing bowl, and mixes them using the pestle. The spatula is used to transfer the material into a commercially available syringe, which is used to inject the material into the bony void. This kit will be provided as a single use, sterile product packaged in a double sterile barrier using industry standard peel pouches.

Skeletal Kinetics' bone void filler is substantially equivalent to other bone void fillers that are commercially available. The following documents were used in making this determination:

Orthopedic Devices Branch ODE. Draft guidance document for the preparation of pre-market notification applications for orthopedic devices—the basic elements. FDA: July 16, 1997.

Restorative Devices Branch ODE. Class II special controls guidance document: resorbable calcium salt bone void filler device: draft guidance for industry and FDA. FDA: February 7, 2002.

OK
unclassified

Skeletal Kinetics, LLC - Confidential


SKS

Please direct any questions concerning this submission to me at 408.366.5002, or Howard Holstein at 202.637.5813. I also request that FDA notify Skeletal Kinetics, LLC of substantial equivalence by sending a facsimile to me at 408.366.1077.

The Company is aware of its obligations to pay user fee upon invoice.

Thank you in advance for your prompt attention to this submission.

Sincerely,



Duran Yetkinler, M.D., Ph.D.

Director of Clinical and Regulatory Affairs

Enclosure

CDRH SUBMISSION COVER SHEET

Date of Submission:**FDA Document Number:**

Section A Type of Submission				
PMA	PMA Supplement	PDP	510(k)	Meeting
<input type="checkbox"/> Original Submission	<input type="checkbox"/> Regular	<input type="checkbox"/> Presubmission Summary	Original Submission	<input type="checkbox"/> Pre-IDE
<input type="checkbox"/> Modular Submission	<input type="checkbox"/> Special	<input type="checkbox"/> Original PDP	<input checked="" type="checkbox"/> Traditional	<input type="checkbox"/> Pre-PMA
<input type="checkbox"/> Amendment	<input type="checkbox"/> Panel Track	<input type="checkbox"/> Notice of intent to start clinical trials	<input type="checkbox"/> Special	<input type="checkbox"/> Pre-PDP
<input type="checkbox"/> Report	<input type="checkbox"/> 30-day Supplement	<input type="checkbox"/> Intention to submit Notice of Completion	<input type="checkbox"/> Abbreviated	<input type="checkbox"/> 180 -day
<input type="checkbox"/> Report Amendment	<input type="checkbox"/> 30-day Notice	<input type="checkbox"/> Notice of Completion	Additional Information	<input type="checkbox"/> Other (specify):
	<input type="checkbox"/> 135-day Supplement	<input type="checkbox"/> Amendment to PDP	<input type="checkbox"/> Traditional	
	<input type="checkbox"/> Real-time Review	<input type="checkbox"/> Report	<input type="checkbox"/> Special	
	<input type="checkbox"/> Amendment to PMA Supplement		<input type="checkbox"/> Abbreviated	
			<input type="checkbox"/> Report Amendment	
IDE	Humanitarian Device Exemption	Class II Exemption	Evaluation of Class III Designation	Other Submission
<input type="checkbox"/> Original Submission	<input type="checkbox"/> Original Submission	<input type="checkbox"/> Original Submission	<input type="checkbox"/> Original Submission	Describe Submission:
<input type="checkbox"/> Amendment	<input type="checkbox"/> Amendment	<input type="checkbox"/> Additional Information	<input type="checkbox"/> Amendment	
<input type="checkbox"/> Supplement	<input type="checkbox"/> Supplement		<input type="checkbox"/> Supplement	
	<input type="checkbox"/> Report			
Section B Applicant or Sponsor				
Company/Institution Name Skeletal Kinetics, LLC		Establishment registration number Skeletal Kinetics has filed an initial device establishment registration form, but FDA has not yet assigned an establishment registration number. Owner/Operator number: 9054231		
Division Name (if applicable):		Phone number (include area code) 408-366-5002		
Street Address 10201 Bubb Road		Fax number (include area code) 408-366-1077		
City Cupertino	State/Province CA	Zip code 95014-4167	Country USA	
Contact Name Duran N. Yetkinler MD, PhD				
Contact Title Director of Clinical and Regulatory Affairs		Contact e-mail address duran@skeletalkinetics.com		
Section C Submission Correspondent (if different from above)				
Company/Institution Name Hogan & Hartson, LLC		Establishment registration number n/a		
Division Name (if applicable):		Phone number (include area code) 202-637-5813		
Street Address 555 13th Street N.W.		Fax number (include area code) 202-637-5910		
City Washington	State/Province DC	Zip code 20004	Country USA	
Contact Name Howard Holstein, Esq.				
Contact Title Regulatory Counsel		Contact e-mail address HMHolstein@hhlaw.com		

Section D1 Reason for Submission --- PMA, PDP or HDE		
<input type="checkbox"/> New Device	<input type="checkbox"/> Change in design, component or specification	<input type="checkbox"/> Location Change
<input type="checkbox"/> Withdrawal	<input type="checkbox"/> Software	<input type="checkbox"/> Manufacturer
<input type="checkbox"/> Additional or expanded indications	<input type="checkbox"/> Color Additive	<input type="checkbox"/> Sterilizer
<input type="checkbox"/> Licensing agreement	<input type="checkbox"/> Material	<input type="checkbox"/> Packager
<input type="checkbox"/> Processing change	<input type="checkbox"/> Specifications	<input type="checkbox"/> Distributor
<input type="checkbox"/> Manufacturing	<input type="checkbox"/> Other (specify below)	<input type="checkbox"/> Report Submission
<input type="checkbox"/> Sterilization	<input type="checkbox"/> Labeling changes	<input type="checkbox"/> Annual or Periodic
<input type="checkbox"/> Packaging	<input type="checkbox"/> Indications	<input type="checkbox"/> Post Approval Study
<input type="checkbox"/> Other (specify below)	<input type="checkbox"/> Instructions	<input type="checkbox"/> Adverse Reaction
<input type="checkbox"/> Response to FDA correspondence	<input type="checkbox"/> Performance Characteristics	<input type="checkbox"/> Device Defect
<input type="checkbox"/> Request for applicant hold	<input type="checkbox"/> Shelf Life	<input type="checkbox"/> Amendment
<input type="checkbox"/> Request for removal of applicant hold	<input type="checkbox"/> Trade Name	<input type="checkbox"/> Change in ownership
<input type="checkbox"/> Request for extension	<input type="checkbox"/> Other (specify below)	<input type="checkbox"/> Change in Correspondent
<input type="checkbox"/> Request to remove or add manufacturing site		
<input type="checkbox"/> Other reason (specify)		

Section D2 Reason for Submission --- IDE		
<input type="checkbox"/> New Device	Change in	<input type="checkbox"/> Response to FDA letter concerning
<input type="checkbox"/> Addition of institution	<input type="checkbox"/> Correspondent	<input type="checkbox"/> Conditional approval
<input type="checkbox"/> Expansion/extension of study	<input type="checkbox"/> Design	<input type="checkbox"/> Deemed approval
<input type="checkbox"/> IRB certification	<input type="checkbox"/> Informed consent	<input type="checkbox"/> Deficient final report
<input type="checkbox"/> Request hearing	<input type="checkbox"/> Manufacturer	<input type="checkbox"/> Deficient progress report
<input type="checkbox"/> Request waiver	<input type="checkbox"/> Manufacturing process	<input type="checkbox"/> Deficient investigator report
<input type="checkbox"/> Termination of study	<input type="checkbox"/> Protocol -- feasibility	<input type="checkbox"/> Disapproval
<input type="checkbox"/> Withdrawal of application	<input type="checkbox"/> Protocol -- other	<input type="checkbox"/> Request extension for time to respond to FDA
<input type="checkbox"/> Unanticipated adverse effect	<input type="checkbox"/> Sponsor	<input type="checkbox"/> Request meeting
<input type="checkbox"/> Notification of emergency use	<input type="checkbox"/> Report Submission	
<input type="checkbox"/> Compassionate use request	<input type="checkbox"/> Current investigator	
<input type="checkbox"/> Treatment IDE	<input type="checkbox"/> Annual progress	
<input type="checkbox"/> Continuing availability request	<input type="checkbox"/> Site waiver limit reached	
<input type="checkbox"/> Other reason (specify)	<input type="checkbox"/> Final	

Section D3 Reason for Submission --- 510(k)		
<input checked="" type="checkbox"/> New Device		
<input type="checkbox"/> Additional or expanded indications	<input type="checkbox"/> Change in Technology	<input type="checkbox"/> Change in materials
<input type="checkbox"/> Other reason (specify)	<input type="checkbox"/> Change in Design	<input type="checkbox"/> Change in manufacturing process

25

Section E Additional Information on 510(k) Submissions				
Product codes of devices to which substantial equivalence is claimed				Summary of, or statement concerning safety and effectiveness data
1 MQV	2	3	4	<input checked="" type="checkbox"/> 510(k) summary attached
5	6	7	8	<input type="checkbox"/> 510(k) statement
510(k) Number	Trade of Proprietary or Model Name			Manufacturer
1 K011897	1 Norian SRS Bone Void Filler			1 Synthes USA
2	2			2
3	3			3
4	4			4
5	5			5
6	6			6

Section F Product Information – Applicable to All Applications					
Common or usual name or classification name Bone Void Filler					
Trade or proprietary or model name			Model Number		
1	Callos™		1	n/a	
2			2		
3			3		
4			4		
5			5		
FDA document numbers of all prior related submissions (regardless of outcome)					
1	2	3	4	5	6
7	8	9	10	11	12
Data included in submission		<input checked="" type="checkbox"/> Laboratory testing	<input checked="" type="checkbox"/> Animal trials	<input type="checkbox"/> Human trials	

Section G Product Classification – Applicable to All Applicants		
Product code: MQV	C.F.R. Section 21 CFR§..... (unclassified)	Device Class
Classification Panel Orthopedic devices branch 87		<input type="checkbox"/> Class I <input type="checkbox"/> Class II <input type="checkbox"/> Class III <input checked="" type="checkbox"/> Unclassified
Indications (from labeling) Callos Bone Void Filler is indicated to fill bony voids or gaps of the skeletal system (i.e. extremities, spine, pelvis). These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. Callos Bone Void Filler is indicated only for bony voids or gaps that are not intrinsic to the stability of the bony structure. The product provides a bone void filler that resorbs and is replaced by bone during the healing process.		

26

Note: submission of this information does not affect the need to submit a 2891 or 2891a Device Establishment Registration form.	FDA Document Number
---	---------------------

Section H Manufacturing/Packaging/Sterilization Sites Relating to a Submission

<input checked="" type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete	FDA Establishment Number	<input checked="" type="checkbox"/> Manufacturer <input type="checkbox"/> Contract manufacturer	<input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager/relabeler
Company/Institution name Skeletal Kinetics, LLC		Establishment registration number Skeletal Kinetics has filed an initial device establishment registration form, but FDA has not yet assigned an establishment registration number Owner/Operator number: 9054231	
Division name (if applicable)		Phone number (include area code) 408-366-5002	
Street address 10201 Bubb Road		FAX number (include area code) 408-366-1077	
City Cupertino	State/Province CA	Zip code 95014-4167	Country USA
Contact name Duran N. Yetkinler MD, PhD			
Contact title Director of Clinical and Regulatory Affairs		Contact e-mail address duran@skeletalkinetics.com	
<input checked="" type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete	FDA Establishment Number (b)	<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract manufacturer	<input checked="" type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager/relabeler

(b) (4)

<input checked="" type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete	FDA Establishment Number (b)(4)	<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract manufacturer	<input checked="" type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager/relabeler
--	------------------------------------	---	--

(b) (4)

27

Table of Contents

I. Administrative Documents..... 8

 1. 510(k) Elements Checklist 8

 2. 510(k) Truthful and Accurate Statement..... 11

 3. 510(k) Indications for Use..... 11

 4. 510(k) Summary..... 11

II. Subject Device..... 12

 1. Device Description..... 12

 a. Device Name 12

 b. Sponsor Registration Number 12

 c. Classification..... 12

 d. Indications for Use 12

 e. Consensus Standards 12

 2. Testing..... 13

 a. Summary of Biocompatibility Testing 13

 b. Summary of Material Characterization Testing 14

 c. Summary of Physical Properties Testing 22

 d. Summary of Performance Testing..... 23

 3. Draft IFU 41

 4. Draft labels 41

 5. Sterilization and shelf life 41

 6. Substantial Equivalence Comparison..... 42

Appendices 45

 Appendix A—510(k) Truthful and Accurate Statement 45

 Appendix B—510(k) Indications for Use 46

 Appendix C—510(k) Summary 47

 Appendix D—Biocompatibility Complete Test Report..... 48

 Appendix E—Material Characterization Complete Test Report 49

 Appendix F—Physical Properties Complete Test Report..... 50

 Appendix G—Performance Complete Test Report 51

 Appendix H—Draft IFU 52

 Appendix I—Draft package labels 53

 Appendix J—Sterilization Validation Protocol..... 54

 Appendix K—LAL Protocol..... 55

28

I. Administrative Documents

1. 510(k) Elements Checklist

Page 1 of 3

510(k) Number: _____

The cover letter clearly identifies the type of 510(k) submission as **(Check the appropriate box)**:

- Special 510(k) - Do Sections 1 and 2
- Abbreviated 510(k) - Do Sections 1, 3 and 4
- Traditional 510(k) or no identification provided - Do Sections 1 and 4

Section 1: Required Elements for All Types of 510(k) submissions:

	Present	Inadequate or Missing
Cover letter, containing the elements listed on page 3-2 of the <u>Premarket Notification [510)] Manual</u> .	p.1	
Table of Contents.	p.7	
Truthful and Accurate Statement.	p.11, (App A)	
Device's Trade Name, Device's Classification Name and Establishment Registration Number.	p.12	
Device Classification Regulation Number and Regulatory Status (Class I, Class II, Class III or Unclassified).	p.12	
Proposed Labeling including the material listed on page 3-4 of the <u>Premarket Notification [510)] Manual</u> .	p.41 (App H & I)	
Statement of Indications for Use that is on a separate page in the premarket submission.	p.11 (App. B)	
Substantial Equivalence Comparison, including comparisons of the new device with the predicate in areas that are listed on page 3-4 of the <u>Premarket Notification [510)] Manual</u> .	p.42	
510(k) Summary or 510(k) Statement.	p.11	

29

	(App. C)	
510(k) Elements Checklist Page 2 of 3		
Description of the device (or modification of the device) including diagrams, engineering drawings, photographs or service manuals.	p.12	
Identification of legally marketed predicate device. *	p.42	
Compliance with performance standards. * [See Section 514 of the Act and 21 CFR 807.87 (d).]	p.12, p.21	
Class III Certification and Summary. **	n/a	
Financial Certification or Disclosure Statement for 510(k) notifications with a clinical study. * [See 21 CFR 807.87 (i)]	n/a	
510(k) Kit Certification ***	n/a	

* - May not be applicable for Special 510(k)s.

** - Required for Class III devices, only.

*** - See pages 3-12 and 3-13 in the Premarket Notification [510] Manual and the Convenience Kits Interim Regulatory Guidance.

Section 2: Required Elements for a SPECIAL 510(k) submission: Not Applicable

Section 3: Required Elements for an ABBREVIATED 510(k)* submission: Not Applicable

Section 4: Additional Requirements for ABBREVIATED and TRADITIONAL 510(k) submissions (If Applicable):

	Present	Inadequate or Missing
a) Biocompatibility data for all patient-contacting materials, OR certification of identical material/formulation:	p.13	
b) Sterilization and expiration dating information:	p.41	
i) sterilization process	p.41	

30

510(k) Elements Checklist Page 3 of 3		
ii) validation method of sterilization process	p.41	
iii) SAL	p.41	
iv) packaging	p.41	
v) specify pyrogen free	p.41	
vi) ETO residues	n/a	
vii) radiation dose	p.41	
c) Software Documentation:	n/a	

Items with checks in the "Present but Deficient" column require additional information from the sponsor. Items with checks in the "Missing" column must be submitted before substantive review of the document.

Passed Screening **Yes** **No**

Reviewer: _____

Concurrence by Review Branch: _____

Date: _____

The deficiencies identified above represent the issues that we believe need to be resolved before our review of your 510(k) submission can be successfully completed. In developing the deficiencies, we carefully considered the statutory criteria as defined in Section 513(i) of the Federal Food, Drug, and Cosmetic Act for determining substantial equivalence of your device. We also considered the burden that may be incurred in your attempt to respond to the deficiencies. We believe that we have considered the least burdensome approach to resolving these issues. If, however, you believe that information is being requested that is not relevant to the regulatory decision or that there is a less burdensome way to resolve the issues, you should follow the procedures outlined in the "A Suggested Approach to Resolving Least Burdensome Issues" document. It is available on our Center web page at: <http://www.fda.gov/cdrh/modact/leastburdensome.html>

31

2. 510(k) Truthful and Accurate Statement

Please see Appendix A for the 510(k) Truthful and Accurate Statement.

3. 510(k) Indications for Use

Please see Appendix B for the 510(k) Indications for Use.

4. 510(k) Summary

Please see Appendix C for the 510(k) Summary

II. Subject Device

1. Device Description

a. Device Name

Callos Bone Void Filler

b. Sponsor Registration Number

Skeletal Kinetics, LLC has applied for an establishment registration number, but FDA has not yet assigned an establishment registration number. However, FDA has informed Skeletal Kinetics, LLC that the Company is considered to be registered. Skeletal Kinetics' Owner/Operator No is 9054231.

c. Classification

FDA has proposed to classify resorbable calcium salt bone void fillers as Class II devices, but the Agency has not finalized that proposed rule. Currently Callos is an unclassified Bone Graft Substitute.

d. Indications for Use

Callos Bone Void Filler is indicated to fill bony voids or gaps of the skeletal system (i.e. extremities, spine, and pelvis). These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. Callos Bone Void Filler is indicated only for bony voids or gaps that are not intrinsic to the stability of the bony structure. The product provides a bone void filler that resorbs and is replaced by bone during the healing process.

e. Consensus Standards

Elemental analysis of Callos Bone Void Filler shows that the device conforms to USP National Formulary (NF) of "Official Monograph of Calcium Sulfate" for applicable items including Iron and Heavy Metal limits, and to USP National Formulary (NF) of "Official Monograph for Tribasic Calcium Phosphate" for applicable items including Arsenic and Heavy Metal limits.

ASTM F1185-88 "Standard Specification for Composition of Ceramic Hydroxyapatite for Surgical Implants" was withdrawn prior to preparation of this 510(k); therefore no comparison was made to this standard.

2. Testing

(b) (4)



34

b. Summary of Material Characterization Testing

(b) (4)



35

(b) (4)



(b) (4)



37

(b) (4)



(b) (4)



(b) (4)



Discussion and Conclusion

(b) (4)



(b) (4)



41

ii. Elemental Analysis

(b) (4)

A large rectangular area of the document is redacted with a solid grey fill, covering the content of the 'Elemental Analysis' section.

Conclusion

(b) (4)

A large rectangular area of the document is redacted with a solid grey fill, covering the content of the 'Conclusion' section.

42

c. Summary of Physical Properties Testing

(b) (4)



i. Physical Characteristics: Density, Porosity, and Dimensional Stability

(b) (4)



ii. Biological Source Material

Callos bone void filler contains no biological source material, e.g., neither animal nor human tissue.

43

d. Summary of Performance Testing

(b) (4)

A large rectangular area of the document is redacted with a solid grey fill, covering the content under the 'Summary of Performance Testing' section.

i. Injection

(b) (4)

A large rectangular area of the document is redacted with a solid grey fill, covering the content under the 'Injection' section.

ii. Setting Time

(b) (4)

A large rectangular area of the document is redacted with a solid grey fill, covering the content under the 'Setting Time' section.

4

(b) (4)



iii. Working time

(b) (4)



45

iv. Temperature and pH measurements

(b) (4)



46

v. ***In vitro* solubility and dissolution testing**

(b) (4)



...phosphate buffered saline.

47

Discussion and Conclusion

(b) (4)



48

vi. Animal testing: An *in vivo* evaluation of two calcium phosphate bone void fillers

(b) (4)



49

(b) (4)



a. Surgical Technique

(b) (4)



(b) (4)



(b) (4)



(b) (4)



d. Chemical and Crystallographic Analysis

(b) (4)



Results

(b) (4)



(b) (4)



53

(b) (4)



(b) (4)



(b) (4)



(b) (4)



b. Biomechanics

(b) (4)



57

Biomechanical Strength Testina

(b) (4)



58

(b) (4)



59

(b) (4)



60

Discussion and Conclusion

(b) (4)



3. Draft IFU

Please see Appendix H for a copy of the draft IFU.

4. Draft labels

Please see Appendix I for sample draft labels for 5 cc Callos Bone Void Filler.

5. Sterilization and shelf life

(b) (4)



(b) (4)



6. Substantial Equivalence Comparison

Table 2: Subject and Predicate Device Comparisons

Substantial Equivalence Comparison	Subject Device Callos Bone Void Filler	Predicate Device Norian SRS Bone Void Filler K011897
Intended Use	(b) (4)	A non-structural bone void filler for osseous defects
Indications for Use		Indicated to fill bony voids or gaps of the skeletal system (i.e. extremities, spine, pelvis). These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. The device is indicated only for bony voids or gaps that are not intrinsic to the stability of the bony structure. The product provides a bone void filler that resorbs and is replaced by bone during the healing process.
Target Population		Individuals with bony defects resulting from surgery or trauma
Design		Self setting calcium phosphate bone void filler which hardens in aqueous environment at 37°C
Components		Package contains one reactant pack with one chamber of powder and one chamber of liquid
Product Preparation		Product is mechanically mixed within a reusable pneumatic mixer
Biocompatibility		Biocompatible
Material Characteristics		
<ul style="list-style-type: none"> • Starting Reactants 		Powder: Alpha tricalcium phosphate, Calcium carbonate (Calcite), Monocalcium phosphate monohydrate. Solution: Dilute sodium phosphate.
<ul style="list-style-type: none"> • Chemical Composition 		Calcium Phosphate Salt
<ul style="list-style-type: none"> • Crystal structure, after hardening 		Hydroxyapatite

63

<ul style="list-style-type: none"> • Calcium to phosphate ratio 	(b) (4)	1.67
Physical Properties		
<ul style="list-style-type: none"> • Porosity 		~50%
<ul style="list-style-type: none"> • Density 		~1.6-1.7 gm/cc
Performance Characteristics		
<ul style="list-style-type: none"> • Injectability 		Injectable for 5 minutes
<ul style="list-style-type: none"> • Working time 		~2 minutes
<ul style="list-style-type: none"> • Setting Time 		~10 minutes
<ul style="list-style-type: none"> • pH 		Physiologic
<ul style="list-style-type: none"> • Setting reaction temperature 		Isothermic
<ul style="list-style-type: none"> • Solubility and Dissolution 		Similar to hydroxyapatite
<ul style="list-style-type: none"> • Bone Remodeling 		New bone grows into the graft area via osteoconduction. The material is replaced by cell mediated remodeling tissue response
Sterility		Sterilized by gamma radiation, single use only
Available Sizes		3cc, 5 cc and 10 cc kits
Voluntary Standards Met	Presumed	

a. Discussion of similarities and differences

Similarities

(b) (4)



64

(b) (4)



Differences and Discussion

(b) (4)



Conclusion

Based on the many similarities and the minimal differences that do not negatively impact the performance of the Callos device, the clinical performance and use of both bone fillers is expected to be essentially indistinguishable. Therefore, Callos bone void filler is substantially equivalent to the predicate device, Norian SRS.

65

Appendices

Appendix A—510(k) Truthful and Accurate Statement

1 page

66

510(k) Truthful and Accurate Statement

[as required by 21 CFR 807.87(j)]

I certify that in my capacity, as the Director of Clinical and Regulatory Affairs for Skeletal Kinetics, LLC, that all data and information submitted in this pre-market notification are truthful and accurate, and that no material facts have been omitted.

Duran Yetkinler

2/20/2003

**SIGNATURE

DATE

Duran Yetkinler, M.D., Ph.D.
Director of Clinical and Regulatory Affairs
Skeletal Kinetics, LLC

PREMARKET NOTIFICATION (510(K)) NUMBER

*As required by 21 CFR Section 807.87(j), effective March 14, 1995.

**Must be signed by a responsible person of the firm required to submit the pre-market notification (e.g., not a consultant for the 510(k) submitter.)

67

Appendix B—510(k) Indications for Use

1 page

68

510(k) Indications for Use

Page ____ of ____

510(k) Number (if known): _____

Device Name: Callos Bone Void Filler

Indications for use:

Callos Bone Void Filler is indicated to fill bony voids or gaps of the skeletal system (i.e. extremities, spine, pelvis). These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. Callos Bone Void Filler is indicated only for bony voids or gaps that are not intrinsic to the stability of the bony structure. The product provides a bone void filler that resorbs and is replaced by bone during the healing process.

(PLEASE DO NOT WRITE BELOW THIS LINE—CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use _____ OR Over-the-Counter Use _____
(per 21 CFR 801.109)

(Division Sign-Off)
Division of General, Restorative and
Neurological Devices

510(k) Number _____

69

Appendix C—510(k) Summary

2 pages

70

510(k) Summary

Sponsor: Skeletal Kinetics, LLC
10201 Bubb Road, Cupertino, CA 95014
Contact Person: Duran Yetkinler, M.D., Ph.D.
Phone Number: 408 366 5002
Fax Number: 408 366 1077
Prepared: February 18, 2003

Trade Name: Callos™
Common Name: Bone Graft Substitute
Classification: Unclassified
Product Code: 87 MVQ

Predicate Device: Callos Bone Void Filler is substantially equivalent to Norian SRS Bone Void Filler (K011897).

Device Description: Callos Bone Void Filler is an injectable, moldable and biocompatible bone void filler. Callos Bone Void Filler resorbs and is replaced with bone during the healing process. The 3 cc, 5 cc, and 10 cc Callos Bone Void Filler kits are provided sterile and are for single use only.

Intended Use/Indications for Use: Callos Bone Void Filler is indicated to fill bony voids or gaps of the skeletal system (i.e. extremities, spine, and pelvis). These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. Callos Bone Void Filler is indicated only for bony voids or gaps that are not intrinsic to the stability of the bony structure. The product provides a bone void filler that resorbs and is replaced by bone during the healing process.

Technological Characteristics: Similar to the predicate device, Callos Bone Void Filler is an injectable, moldable, biocompatible, resorbable calcium phosphate based material intended for identical indications.

Performance Data: Non-clinical testing included material properties such as density, porosity, dimensional stability, injectability, setting time, working time, pH, and setting temperature. Biocompatibility testing demonstrated that the material is non-cytotoxic, non-systemic toxic, non-mutagenic, non-irritative, non-pyrogenic, and non-sensitizing. Comparative testing with the predicate device showed equivalence in terms of solubility and dissolution rate, X-Ray Diffraction (XRD), Fourier Transform Infrared (FTIR) spectroscopy and elemental analysis. Animal testing demonstrated substantial equivalence to the predicate device following *in vivo* implantation. Histological, chemical, crystallographical, and mechanical analyses showed substantial equivalence.

Basis for Substantial Equivalence: The Callos Bone Void Filler has the same intended use, identical indications, and very similar technological characteristics as the predicate device. Any minor technological differences between Callos Bone Void Filler and its predicate device do not raise any new issues of safety or effectiveness.

Functional, biocompatibility, and animal testing results show that the Callos Bone Void Filler is as safe and effective as the predicate device. Thus, the Callos Bone Void Filler is substantially equivalent.

72

Appendix D—Biocompatibility Complete Test Report

1. MEM Elution Test: <i>In vitro</i> cytotoxicity:	8 pages
2. <i>Salmonella Typhimurium</i> Reverse Mutation Assay: AMES Test	17 pages
3. Irritation: Intracutaneous reactivity (ISO)	9 pages
4. Systemic Toxicity: ISO/USP Systemic Injection	8 pages
5. Sensitization: ISO Magnusson Klingman Method (2 extracts)	13 pages

73

Appendix E—Material Characterization Complete Test Report

1. Chemical and Crystallographic analysis: Fourier transform infrared (FTIR) and X-ray diffraction (XRD) spectroscopy 8 pages
2. Elemental Analysis 6 pages

Title: **Chemical and Crystallographic analysis:** (b) (4)

Investigators: Dave Delaney, Brent R Constantz PhD

Approval: Duran N Yetkinler MD, PhD

Facility: Skeletal Kinetics, LLC
10211 Bubb Road
Cupertino, CA 95014-4166

Testing: 18-January-03

Final Report: 31-January-03

Objective

(b) (4)



(b) (4)



Results

(b) (4)



131

(b) (4)



132

(b) (4)



133

Results

(b) (4)



134

(b) (4)



138

(b) (4)



Discussion and Conclusion

(b) (4)



36

(b) (4)



137

Appendix F—Physical Properties Complete Test Report

Physical Characteristics: Density, Porosity, and Dimensional Stability 2 pages

144

Title: Physical Characteristics: Density, Porosity, and Dimensional Stability

Investigator: David Delaney

Approval: Duran N Yetkinler MD, PhD

Facility: Skeletal Kinetics
10201 Bubb Rd
Cupertino CA 95014

Objective (b) (4)

[Redacted]

Materials

(b) (4)

[Redacted]

Density and Porosity Procedure

(b) (4)

[Redacted]

Dimensional Stability Procedure

(b) (4)

[Redacted]

Results

(b) (4)

[Redacted]

145'

(b) (4) [Redacted]

[Redacted]

(b) (4) [Redacted]

(b) (4) [Redacted]

Discussion and Conclusions

(b) (4) [Redacted]

146

Appendix G—Performance Complete Test Report

1. Injection Testing of Callos Bone Void Filler	3 pages
2. Initiation of Setting Testing of Callos Bone Void Filler	3 pages
3. Working Time Testing of Callos Bone Void Filler	3 pages
4. Temperature and pH Testing of Callos Bone Void Filler	3 pages
5. <i>In vitro</i> Solubility and Dissolution Rate	4 pages
6. Animal Study: An <i>in vivo</i> evaluation of two calcium phosphate bone void fillers	14 pages

147

Title: Injection Testing of Callos Bone Void Filler

Investigators: Dave Delaney, Brent R Constantz PhD

Approval: Duran N Yetkinler MD, PhD

Facility: Skeletal Kinetics, LLC
10211 Bubb Road
Cupertino, CA 95014-4166

Testing: 18-January-03

Final Report: 31-January-03

(b) (4)



(b) (4)



Results

(b) (4)



Discussion and Conclusion

Skeletal Kinetics, LLC - Confidential

Page 2 of 3 Appendix G.1

149

(b) (4)



Table 1: Injection Test Data

(b) (4)



References

ASTM F 451-86

150

Title: **Initiation of Setting Testing of Callos Bone Void Filler**

Investigators: Dave Delaney, Brent R Constantz PhD

Approval: Duran N Yetkinler MD PhD

Facility: Skeletal Kinetics, LLC
 10211 Bubb Road
 Cupertino, CA 95014

Testing: 22-January-03

Final Report: 31-January-03

(b) (4)



(b) (4)



Results

(b) (4)



Conclusions

(b) (4)



References

ASTM C403/C403M-99
ASTM C266-99

152

Table 1: (b) (4)

(b) (4)



153

Title: Working Time Testing of Callos Bone Void Filler

Investigators: David Delaney, Brent R Constantz PhD

Approval: Duran N Yetkinler MD PhD

Facility: Skeletal Kinetics, LLC
10211 Bubb Road
Cupertino, CA 95014

Testing: 20-January-03

Final Report: 31-January-03

(b) (4)



Procedure

(b) (4)



Results

(b) (4)



Conclusions

(b) (4)



References

ASTM C403/C403M-99
ASTM C266-99

155-

(b) (4)



156

Title: Temperature and pH Testing of Callos Bone Void Filler

Investigator: David Delaney

Approval: Duran N Yetkinler MD, PhD

Facility: Skeletal Kinetics
10201 Bubb Rd
Cupertino CA 95014

(b) (4)



Results

(b) (4)



(b) (4)



158

Discussion and Conclusion

(b) (4)



159

Title: *In vitro* Solubility and Dissolution Rate

Investigators: Dave Delaney, Brent R Constantz PhD

Approval: Duran N Yetkinler MD, PhD

Facility: Skeletal Kinetics, LLC
10211 Bubb Road
Cupertino, CA 95014-4166

Testing: 18-January-03

Final Report: 31-January-03

(b) (4)



(b) (4)

A large rectangular area of the document is completely redacted with a solid grey fill.

Results

(b) (4)

A large rectangular area of the document is completely redacted with a solid grey fill, covering the entire 'Results' section.

161

(b) (4)



Discussion and Conclusions

(b) (4)



162

(b) (4)



163

AN *IN VIVO* EVALUATION OF TWO CALCIUM PHOSPHATE BONE VOID FILLERS

PURPOSE

This study evaluated the biocompatibility, bone resorption, formation, and biomechanical properties of the implanted regions containing Callos and Norian SRS bone void fillers. Histological, mechanical and crystallographic analyses were performed following implantation in the rabbit femur. Callos™ (Skeletal Kinetics, LLC, Cupertino, CA) and Norian SRS (Synthes Corp., Paoli, PA) were used in this paired design study.

METHODS

(b) (4)



164

Table 1. Flow chart of the experimental design for the study.

(b) (4)



Surgical Technique

(b) (4)



165

(b) (4)



Histological Analysis

(b) (4)



166

(b) (4)



Biomechanical Testing

(b) (4)



Chemical and Crystallographic Analysis

Skeletal Kinetics, LLC - Confidential

Page 4 of 14

Appendix G.6

167

(b) (4)



Statistical Analysis

(b) (4)



RESULTS

Histology

(b) (4)



(b) (4)



168

(b) (4)



(b) (4)



170

(b) (4)



171

(b) (4)



172

Biomechanics

(b) (4)



173

Chemistry and Crystallography

(b) (4)



174

(b) (4)



175

DISCUSSION AND CONCLUSION

(b) (4)



176

REFERENCES

(b) (4)



177

Appendix H—Draft IFU

3 pages

Draft IFU

(b) (4)



9

(b) (4)



10. Instructions for Use:

(b) (4)



181

Appendix I—Draft package labels

2 pages

182

Draft package labels

(b) (4)



(b) (4)



183

(b) (4)



184

Appendix J—Sterilization Validation Protocol

8 pages

185

DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service
Food and Drug Administration

Memorandum

From: Reviewer(s) - Name(s)

William Sloan

Subject: 510(k) Number

K030554

To: The Record - It is my recommendation that the subject 510(k) Notification:

Refused to accept.

Requires additional information (other than refuse to accept).

Is substantially equivalent to marketed devices.

NOT substantially equivalent to marketed devices.

De Novo Classification Candidate?

YES

NO

Other (e.g., exempt by regulation, not a device, duplicate, etc.)

Is this device subject to Postmarket Surveillance?

YES

NO

Is this device subject to the Tracking Regulation?

YES

NO

Was clinical data necessary to support the review of this 510(k)?

YES

NO

Is this a prescription device?

YES

NO

Was this 510(k) reviewed by a Third Party?

YES

NO

Special 510(k)?

YES

NO

Abbreviated 510(k)? Please fill out form on H Drive 510k/boilers

YES

NO

This 510(k) contains:

Truthful and Accurate Statement Requested Enclosed
(required for originals received 3-14-95 and after)

A 510(k) summary OR A 510(k) statement

The required certification and summary for class III devices

The indication for use form (required for originals received 1-1-96 and after)

Animal Tissue Source YES NO

COMBINATION: N

The submitter requests under 21 CFR 807.95 (doesn't apply for SEs):

No Confidentiality Confidentiality for 90 days Continued Confidentiality exceeding 90 day

Predicate Product Code with class:

Additional Product Code(s) with panel (optional):

now unclassified

Review:
(Branch Chief)

REDB
(Branch Code)

5/19/03
(Date)

Final Review:
(Division Director)

Mark N. Millman


5/19/03
(Date)

Revised: 8/17/99

fr

4

510(k) "SUBSTANTIAL EQUIVALENCE" DECISION-MAKING PROCESS



- ❖ 510(k) Submissions compare new devices to marketed devices. FDA requests additional information if the relationship between marketed and "predicate" (pre-Amendments or reclassified post-Amendments) devices is unclear.
- ❖❖ This decision is normally based on descriptive information alone, but limited testing information is sometimes required.
- ❖❖❖ Data maybe in the 510(k), other 510(k)s, the Center's classification files, or the literature.

5

"SUBSTANTIAL EQUIVALENCE" (SE) DECISION MAKING DOCUMENTATION

K030554

Date: 5/16/03

Reviewer: Nadine Y. Sloan

Division/Branch: DGRND/REDB

Device Name: Callos Bone Void Filler

Product To Which Compared (510(K) Number If Known: Norian SRS (K011897)

		YES	NO	
1.	Is Product A Device	✓		If NO = Stop
2.	Is Device Subject To 510(k)?	✓		If NO = Stop
3.	Same Indication Statement?	✓		If YES = Go To 5
4.	Do Differences Alter The Effect Or Raise New Issues Safety Or Effectiveness?			If YES = Stop NE
5.	Same Technological Characteristics?		✓	If YES = Go To 7
6.	Could The New Characteristics Affect Safety Or Effectiveness?		✓	If YES = Go To 8
7.	Descriptive Characteristics Precise Enough?		✓	If NO = Go To 10 If YES = Stop SE
8.	New Types Of Safety Or Effectiveness Questions?			If YES = Stop NE
9.	Accepted Scientific Methods Exist?			If NO = Stop NE
10.	Performance Data Available?	✓		If NO = Request Data
11.	Data Demonstrate Equivalence?	✓		Final Decision: SE

Note: In addition to completing the form on the LAN, "yes" responses to questions 4, 6, 8, and 11, and every "no" response requires an explanation.

1. **Intended Use:** Identical to predicate device.

Callos Bone Void filler is indicated to fill bony voids or gaps of the skeletal system (i.e., extremities, spine and pelvis). These defects may be surgically created osseous defects or osseous defects created from traumatic injury to bone. Callos Bone Void filler is indicated only for bony voids or gaps that are not intrinsic to the stability of the bony structure. The product provides a bone void filler that resorbs and is replaced with bone during the healing process.

2. **Device Description:** Provide a statement of how the device is either similar to and/or different from other marketed devices, plus data (if necessary) to support the statement.

Subject Device: (b) (4)

[Redacted]

[Redacted]

[Redacted]

Predicate Device: *Norian is a calcium phosphate based bone cement which hardens in situ to form a carbonated apatite. The pre-mixed Norian SRS Bone Void Filler is comprised of monocalcium phosphate, monohydrate [MCPM, (Ca(H₂PO₄)₂.H₂O)], α-tricalcium phosphate [TCP, Ca₃(PO₄)₂], and calcium carbonate (CC, CaCO₃), to which a sodium phosphate solution (NaHPO₄.7H₂O) is added to form a paste. The device is injected into the void space where it hardens in approximately 10 minutes. Norian SRS Bone Void Filler is slowly resorbed/remodeled. See attached summary table provided by the sponsor.*

Is the device life-supporting or life sustaining? *No.*

Is the device implanted (short-term or long-term)? *Yes (long term implant that slowly resorbs over a period of years)*

Does the device design use software? *NA*

Is the device sterile? *Yes.*

Is the device for single use? *Yes.*

Is the device over-the-counter or prescription use? *Prescription use.*

Does the device contain drug or biological product as a component? *No.*

Is this device a kit? *No (not in the regulatory sense, although it is provided with components necessary to mix materials to form intended paste).*

Provide a summary about the devices design, materials, physical properties and toxicology profile if important.

(b) (4)

[Redacted]

(b) (4)



comparison table provided by sponsor for results).

EXPLANATIONS TO "YES" AND "NO" ANSWERS TO QUESTIONS ON PAGE 1 AS NEEDED

1. Explain why not a device: NA
2. Explain why not subject to 510(k): NA
3. How does the new indication differ from the predicate device's indication: NA
4. Explain why there is or is not a new effect or safety or effectiveness issue: NA
5. **Describe the new technological characteristics:**

Different starting materials are used, as specified above.

6. **Explain how new characteristics could or could not affect safety or effectiveness:**

(b) (4)



7. **Explain how descriptive characteristics are not precise enough:**

(b) (4)



11. **Explain how the performance data demonstrates that the device is or is not substantially equivalent:**

(b) (4)



ATTACH ADDITIONAL SUPPORTING INFORMATION


See attached SE comparison table and discussion (taken directly from the 510(k)).

INSTRUCTIONS FOR USE

The labeling has been revised to be consistent with the cautions and precautions section for the predicate device. (Not all of the cautions/precautions for the predicate device were incorporated into the Callos labeling since not all of the cautions applied to devices of this type.) The final draft labeling was reviewed in comparison to predicate device and was determined to be acceptable. In addition, the instructions are consistent with how the device was studied (in vitro and in vivo).

RECOMMENDATION:

Substantially equivalent to Norian Bone Void Filler (K011897), unclassified.


Nadine Y. Sloan
Callos.SE.doc

6. Substantial Equivalence Comparison

Table 2: Subject and Predicate Device Comparisons

Substantial Equivalence Comparison	Subject Device Callos Bone Void Filler	Predicate Device Norian SRS Bone Void Filler K011897
Intended Use	(b) (4)	A non-structural bone void filler for osseous defects
Indications for Use		Indicated to fill bony voids or gaps of the skeletal system (i.e. extremities, spine, pelvis). These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. The device is indicated only for bony voids or gaps that are not intrinsic to the stability of the bony structure. The product provides a bone void filler that resorbs and is replaced by bone during the healing process.
Target Population		Individuals with bony defects resulting from surgery or trauma
Design		Self setting calcium phosphate bone void filler which hardens in aqueous environment at 37°C
Components		Package contains one reactant pack with one chamber of powder and one chamber of liquid
Product Preparation		Product is mechanically mixed within a reusable pneumatic mixer
Biocompatibility		Biocompatible
Material Characteristics		
<ul style="list-style-type: none"> <li data-bbox="375 1682 634 1713">• Starting Reactants 		Powder: Alpha tricalcium phosphate, Calcium carbonate (Calcite), Monocalcium phosphate monohydrate. Solution: Dilute sodium phosphate.
<ul style="list-style-type: none"> <li data-bbox="375 1808 566 1866">• Chemical Composition 		Calcium Phosphate Salt
<ul style="list-style-type: none"> <li data-bbox="375 1871 626 1927">• Crystal structure, after hardening 		Hydroxyapatite

10

<ul style="list-style-type: none"> • Calcium to phosphate ratio 	(b) (4)	1.67
Physical Properties		
<ul style="list-style-type: none"> • Porosity 		~50%
<ul style="list-style-type: none"> • Density 		~1.6-1.7 gm/cc
Performance Characteristics		
<ul style="list-style-type: none"> • Injectability 		Injectable for 5 minutes
<ul style="list-style-type: none"> • Working time 		~2 minutes
<ul style="list-style-type: none"> • Setting Time 		~10 minutes
<ul style="list-style-type: none"> • pH 		Physiologic
<ul style="list-style-type: none"> • Setting reaction temperature 		Isothermic
<ul style="list-style-type: none"> • Solubility and Dissolution 		Similar to hydroxyapatite
<ul style="list-style-type: none"> • Bone Remodeling 		New bone grows into the graft area via osteoconduction. The material is replaced by cell mediated remodeling tissue response
Sterility		Sterilized by gamma radiation, single use only
Available Sizes		3cc, 5 cc and 10 cc kits
Voluntary Standards Met	Presumed	

a. Discussion of similarities and differences

Similarities

(b) (4)

11

(b) (4)



Differences and Discussion

(b) (4)



Conclusion

Based on the many similarities and the minimal differences that do not negatively impact the performance of the Callos device, the clinical performance and use of both bone fillers is expected to be essentially indistinguishable. Therefore, Callos bone void filler is substantially equivalent to the predicate device, Norian SRS.

12

Sloan, Nadine Y.

From: duran yetkinler [duran@skeletalkinetics.com]
Sent: Wednesday, May 14, 2003 5:54 PM
To: NYR@CDRH.FDA.GOV
Subject: Changes in the Appendix H

Nadine,

It was nice speaking with you on the phone. I am forwarding the revised draft IFU to reflect the changes that we talked on the phone. I made changes to:

(b) (4)



Let me know if there are any additional questions.

Regards,

Duran

Duran N Yetkinler MD PhD
VP of Regulatory and Product Development
Skeletal Kinetics, LLC
10201 Bubb Road
Cupertino, CA 95014
Phone 408 366 5002
Fax 408 366 1077
Cell 408 757 6603

13

5/15/2003

Draft IFU

(b) (4)



(b) (4)



15

(b) (4)

A large rectangular area of the document is redacted with a solid grey fill, covering the majority of the upper half of the page.

10. Instructions for Use:

(b) (4)

A very large rectangular area of the document is redacted with a solid grey fill, covering the entire lower half of the page and extending up to the 'Instructions for Use' section header.

16

**SCREENING CHECKLIST
FOR ALL PREMARKET NOTIFICATION [510(k)] SUBMISSIONS.**

510(k) Number: K 030554

The cover letter clearly identifies the type of 510(k) submission as (Check the appropriate box):

- Special 510(k) - Do Sections 1 and 2
- Abbreviated 510(k) - Do Sections 1, 3 and 4
- Traditional 510(k) or no identification provided - Do Sections 1 and 4

Section 1: Required Elements for All Types of 510(k) submissions:

	Present or Adequate	Missing or Inadequate
Cover letter, containing the elements listed on page 3-2 of the Premarket Notification [510] Manual.		
Table of Contents.		
Truthful and Accurate Statement.		
Device's Trade Name, Device's Classification Name and Establishment Registration Number.		
Device Classification Regulation Number and Regulatory Status (Class I, Class II, Class III or Unclassified).		
Proposed Labeling including the material listed on page 3-4 of the Premarket Notification [510] Manual.		
Statement of Indications for Use that is on a separate page in the premarket submission.		
Substantial Equivalence Comparison, including comparisons of the new device with the predicate in areas that are listed on page 3-4 of the Premarket Notification [510] Manual.		
510(k) Summary or 510(k) Statement.		
Description of the device (or modification of the device) including diagrams, engineering drawings, photographs or service manuals.		
Identification of legally marketed predicate device. *		
Compliance with performance standards. * [See Section 514 of the Act and 21 CFR 807.87 (d).]		
Class III Certification and Summary. **		
Financial Certification or Disclosure Statement for 510(k) notifications with a clinical study. * [See 21 CFR 807.87 (i)]		
510(k) Kit Certification ***		

not present
 see or copy manuals
 completely for them

* - May not be applicable for Special 510(k)s.
 ** - Required for Class III devices, only.
 *** - See pages 3-12 and 3-13 in the Premarket Notification [510] Manual and the Convenience Kits Interim Regulatory Guidance.

Section 2: Required Elements for a SPECIAL 510(k) submission:

	Present	Inadequate or Missing
Name and 510(k) number of the submitter's own, unmodified predicate device.		
A description of the modified device and a comparison to the sponsor's predicate device.		
A statement that the intended use(s) and indications of the modified device, as described in its labeling are the same as the intended uses and indications for the submitter's unmodified predicate device.		
Reviewer's confirmation that the modification has not altered the fundamental scientific technology of the submitter's predicate device.		
A Design Control Activities Summary that includes the following elements (a-c):		
a. Identification of Risk Analysis method(s) used to assess the impact of the modification on the device and its components, and the results of the analysis.		
b. Based on the Risk Analysis, an identification of the required verification and validation activities, including the methods or tests used and the acceptance criteria to be applied.		
c. A Declaration of Conformity with design controls that includes the following statements:		
A statement that, as required by the risk analysis, all verification and validation activities were performed by the designated individual(s) and the results of the activities demonstrated that the predetermined acceptance criteria were met. This statement is signed by the individual responsible for those particular activities.		
A statement that the manufacturing facility is in conformance with the design control procedure requirements as specified in 21 CFR 820.30 and the records are available for review. This statement is signed by the individual responsible for those particular activities.		

Section 3: Required Elements for an ABBREVIATED 510(k)* submission:

	Present	Inadequate or Missing
For a submission, which relies on a guidance document and/or special control(s), a summary report that describes how the guidance and/or special control(s) was used to address the risks associated with the particular device type. (If a manufacturer elects to use an alternate approach to address a particular risk, sufficient detail should be provided to justify that approach.)		
For a submission, which relies on a recognized standard, a declaration of conformity [For a listing of the required elements of a declaration of conformity, SEE Required Elements for a Declaration of Conformity to a Recognized Standard, which		

is posted with the 510(k) boilers on the H drive.)		
For a submission, which relies on a recognized standard without a declaration of conformity, a statement that the manufacturer intends to conform to a recognized standard and that supporting data will be available before marketing the device.		
For a submission, which relies on a non-recognized standard that has been historically accepted by FDA, a statement that the manufacturer intends to conform to a recognized standard and that supporting data will be available before marketing the device.		
For a submission, which relies on a non-recognized standard that has <u>not</u> been historically accepted by FDA, a statement that the manufacturer intends to conform to a recognized standard and that supporting data will be available before marketing the device <u>and</u> any additional information requested by the reviewer in order to determine substantial equivalence.		
Any additional information, which is not covered by the guidance document, special control, recognized standard and/or non-recognized standard, in order to determine substantial equivalence.		

- * - When completing the review of an abbreviated 510(k), please fill out an Abbreviated Standards Data Form (located on the H drive) and list all the guidance documents, special controls, recognized standards and/or non-recognized standards, which were noted by the sponsor.

Section 4: Additional Requirements for ABBREVIATED and TRADITIONAL 510(k) submissions (If Applicable):

	Present	Inadequate or Missing
a) Biocompatibility data for all patient-contacting materials, OR certification of identical material/formulation:	/	
b) Sterilization and expiration dating information:	/	
i) sterilization process	/	
ii) validation method of sterilization process	/	
iii) SAL	/	
iv) packaging	/	
v) specify pyrogen free	/	
vi) ETO residues	Ag	
vii) radiation dose	/	
viii) Traditional Method or Non-Traditional Method	/	
c) Software Documentation:	n/a	

Items with checks in the "Present or Adequate" column do not require e additional information from the sponsor. Items with checks in the "Missing or Inadequate" column must be submitted before substantive review of the document.

Passed Screening Yes No
 Reviewer: Miriam Swan
 Concurrence by Review Branch: _____

Internal Administrative Form

	YES	NO
1. Did the firm request expedited review?		✓
2. Did we grant expedited review?		
3. Have you verified that the Document is labeled Class III for GMP purposes?	NA	
4. If, not, has POS been notified?		
5. Is the product a device?	✓	
6. Is the device exempt from 510(k) by regulation or policy?	✓	✓
7. Is the device subject to review by CDRH?	✓	
8. Are you aware that this device has been the subject of a previous NSE decision?		✓
9. If yes, does this new 510(k) address the NSE issue(s), (e.g., performance data)?		
10. Are you aware of the submitter being the subject of an integrity investigation?		✓
11. If, yes, consult the ODE Integrity Officer.		
12. Has the ODE Integrity Officer given permission to proceed with the review? (Blue Book Memo #191-2 and Federal Register 90N0332, September 10, 1991.		