510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

Cleo™ 90 Infusion Set

August 10, 2004

OCT 7 - 2004

GENERAL INFORMATION

Applicant's Name and Address: Smiths Medical MD, Inc.
1265 Grey Fox Road
St. Paul, MN 55112

Contact Person: David H. Short
Director Regulatory and Clinical Affairs

Common/Usual Name: Subcutaneous Infusion Set

Proprietary Name: Cleo™ 90 Infusion Set

Equivalence Device Comparison: Unomedical Comfort Set; Disetronic® Ultraflex™ Set; and the Medtronic MiniMed Sof-Set®

II. DEVICE DESCRIPTION

The Cleo™ 90 Infusion Set consists of an applicator that contains an adhesive site with cannula and needle and separate extension tubing set with buckle. The infusion set has a standard luer connection. After insertion of the needle and cannula into the subcutaneous tissue, the applicator automatically retracts the needle into the applicator housing providing needle protection. The infusion set buckle attaches to the site allowing infusion of medication through the extension tubing.

III. INTENDED USE OF THE DEVICE

The Cleo 90 Infusion Set is designed for subcutaneous infusion of medication (including insulin) only. The set is not intended for intravenous (IV) infusion or to infuse blood or blood products.

IV. DEVICE COMPARISON

The Cleo 90 Infusion Set was compared to the following similar products: Unomedical, Comfort Set; Disetronic®, Ultraflex™ Set; and the Medtronic MiniMed, Sof-Set®. The Cleo 90 Infusion Set is substantially equivalent to these products by having the identical indication for use, a site with flexible catheter and needle for insertion into the subcutaneous tissue and a separate extension tubing with buckle for attachment to the site. The Unomedical Comfort Set utilizes a rigid needle rather than a flexible cannula. All sets have an adhesive patch that secures the catheter or needle to the skin following placement. The

042172

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

Cleo™ 90 Infusion Set

August 10, 2004

Unomedical Comfort Set and Medtronic MiniMed Sof-Set have detachable extension sets, where as the Disetronic, Ultraflex Set has an integral extension set. For those devices with a flexible cannula, the needle is removed prior to infusion. The Cleo 90 Infusion Set incorporates a needle protection feature that retracts the needle automatically following insertion of the needle and cannula. The needle may also be retracted automatically by the user as well.

V. SUMMARY OF STUDIES

A. Functional Testing

In-vitro functional testing of the Cleo 90 Infusion Set was conducted. Biocompatibility testing was performed on the materials of the Cleo 90 Infusion Set.

B. Clinical Studies

Human clinical studies were not deemed necessary to evaluate the safety or effectiveness of the Cleo 90 Infusion Set.

To evaluate the safety and effectiveness of the needle protection feature of the Cleo 90 Infusion Set, a simulated clinical use evaluation was conducted in accordance with the "Supplementary Guidance on Premarket Notifications for Medical Devices with Sharps Injury Prevention Features; Guidance for Industry and FDA", December 31, 2002. Individuals who participated in the evaluation included patients who routinely use an insulin pump and infusion set, parents or guardians of patients who insert an infusion set on behalf of a patient, and health care providers who instruct patients and caregivers on the use of infusion pumps and infusion sets. This evaluation determined the needle protection feature operated as intended and no needle sticks or injuries were reported in the evaluation by the participants.

C. Conclusions Drawn from the Studies

The results of the testing conducted indicate the Cleo 90 Infusion Set functioned according to specifications and the materials used in the device are biocompatible. Based on these results, the product is considered acceptable for human use.

DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

OCT 7 - 2004

Mr. David H. Short
Director Regulatory and Clinical Affairs
Smiths Medical MD, Incorporated
1265 Grey Fox Road
St. Paul, Minnesota 55112

Re: K042172
Trade/Device Name: Cleo 90 Infusion Set
Regulation Number: 880.5440
Regulation Name: Intravascular Administration Set
Regulatory Class: II
Product Code: FPA
Dated: August 10, 2004
Received: August 11, 2004

Dear Mr. Short:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Page 2 – Mr. Short

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

K042172

Indications for Use

510(k) Number (if known): K042172

Device Name: Cleo 90 Infusion Set

Indications For Use:

The Cleo 90 Infusion Set is designed for subcutaneous infusion of medication (including insulin) only. The set is not intended for intravenous (IV) infusion or to infuse blood or blood products

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Neola Hubbard for Anthony Watson

(Division Sign-Off)
Division of Anesthesiology, General Hospital,
Infection Control, Dental Devices

510(k) Number: K042172

DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

OCT 7 - 2004

Mr. David H. Short
Director Regulatory and Clinical Affairs
Smiths Medical MD, Incorporated
1265 Grey Fox Road
St. Paul, Minnesota 55112

Re: K042172
Trade/Device Name: Cleo 90 Infusion Set
Regulation Number: 880.5440
Regulation Name: Intravascular Administration Set
Regulatory Class: II
Product Code: FPA
Dated: August 10, 2004
Received: August 11, 2004

Dear Mr. Short:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Page 2 -- Mr. Short

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

K042172

Indications for Use

510(k) Number (if known): K042172

Device Name: Cleo 90 Infusion Set

Indications For Use:

The Cleo 90 Infusion Set is designed for subcutaneous infusion of medication (including insulin) only. The set is not intended for intravenous (IV) infusion or to infuse blood or blood products

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Neola Hubbard for Anthony Watson

(Division Sign-Off)
Division of Anesthesiology, General Hospital,
Infection Control, Dental Devices

510(k) Number: K042172

Page 1 of 1

DEPARTMENT OF HEALTH AND HUMAN SERVICES

Public Health Service

Food and Drug Administration
 Center for Devices and
 Radiological Health
 Office of Device Evaluation
 Document Mail Center (HFZ-401)
 9200 Corporate Blvd.
 Rockville, Maryland 20850

August 11, 2004

SMITHS MEDICAL MD, INC.
 1265 GREY FOX ROAD
 ST. PAUL, MN 55112
 ATTN: DAVID H. SHORT

510(k) Number: K042172
 Received: 11-AUG-2004
 Product: CLEO 90 INFUSION SET

The Food and Drug Administration (FDA), Center for Devices and Radiological Health (CDRH), has received the Premarket Notification you submitted in accordance with Section 510(k) of the Federal Food, Drug, and Cosmetic Act (Act) for the above referenced product. We have assigned your submission a unique 510(k) number that is cited above. Please refer prominently to this 510(k) number in any future correspondence that relates to this submission. We will notify you when the processing of your premarket notification has been completed or if any additional information is required. YOU MAY NOT PLACE THIS DEVICE INTO COMMERCIAL DISTRIBUTION UNTIL YOU RECEIVE A LETTER FROM FDA ALLOWING YOU TO DO SO.

On May 21, 2004, FDA issued a Guidance for Industry and FDA Staff entitled, "FDA and Industry Actions on Premarket Notification (510(k)) Submissions: Effect on FDA Review Clock and Performance Assessment". The purpose of this document is to assist agency staff and the device industry in understanding how various FDA and industry actions that may be taken on 510(k)s should affect the review clock for purposes of meeting the Medical Device User Fee and Modernization Act. Please review this document at <http://www.fda.gov/cdrh/mdufma/guidance/1219.html>.

Please remember that all correspondence concerning your submission MUST be sent to the Document Mail Center (DMC)(HFZ-401) at the above letterhead address. Correspondence sent to any address other than the one above will not be considered as part of your official premarket notification submission. Also, please note the new Blue Book Memorandum regarding Fax and E-mail Policy entitled, "Fax and E-Mail Communication with Industry about Premarket Files Under Review". Please refer to this guidance for information on current fax and e-mail practices at www.fda.gov/cdrh/ode/a02-01.html.

You should be familiar with the regulatory requirements for medical device available at Device Advice <http://www.fda.gov/cdrh/devadvice/>". If you have other procedural or policy questions, or want information on how to check on the status of your submission, please contact DSMICA at (301) 443-6597 or its toll-free number (800) 638-2041, or at their Internet address <http://www.fda.gov/cdrh/dsmamain.html> or me at (301)594-1190.

Sincerely yours,

Marjorie Shulman
 Supervisory Consumer Safety Officer
 Office of Device Evaluation
 Center for Devices and Radiological Health

K042172

smiths

Smiths Medical MD, Inc.

**Medication Delivery
and Patient Monitoring Division**

1265 Grey Fox Road
St Paul MN 55112 USA
Tel: +1 651-633-2556 Fax: +1 651-628-7459
www.deltec.com

August 10, 2004

Food and Drug Administration
Office of Device Evaluation
Document Mail Center (HFZ-401)
Center for Devices and Radiological Health
9200 Corporate Blvd.
Rockville, MD 20850

RE: 510(k) Notification: Cleo™ 90 Infusion Set

Dear Sir or Madam:

Enclosed is a 510(k) Notification for the Cleo 90™ Infusion Set manufactured by Smiths Medical MD, Inc.

Please find attached CDRH Premarket Review Cover Sheet. In addition, a Medical Device User Fee Cover Sheet is also attached:

(PAYMENT IDENTIFICATION NUMBER (b) (4))

This submission is submitted in accordance with 21 CFR, Part 807, Subpart E and is supplied in duplicate as required by regulation.

Sincerely,
Smiths Medical, Inc.

David H. Short
Director Regulatory and Clinical Affairs

Enclosure

HC II

51214

Form Approved OMB No. 0910-0511 Expiration Date August 31, 2006 See instructions for OMB Statement.

DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION MEDICAL DEVICE USER FEE COVER SHEET		PAYMENT IDENTIFICATION NUMBER: (b) (4) Write the Payment Identification Number on your check.
A completed Cover Sheet must accompany each original application or supplement subject to fees. The following actions must be taken to properly submit your application and fee payment:		
<ol style="list-style-type: none"> 1. Electronically submit the completed Cover Sheet to the Food and Drug Administration (FDA) before payment is sent. 2. Include a printed copy of this completed Cover Sheet with a check made payable to the Food and Drug Administration. Remember that the Payment Identification Number must be written on the check. 3. Mail Check and Cover Sheet to the US Bank Lock Box, FDA Account, P.O. Box 956733, St. Louis, MO 63195-6733. (Note: In no case should payment be submitted with the application.) 4. If you prefer to send a check by a courier, the courier may deliver the check and Cover Sheet to: US Bank, Attn: Government Lockbox 956733, 1005 Convention Plaza, St. Louis, MO 63101. (Note: This address is for courier delivery only. Contact the US Bank at 314-418-4821 if you have any questions concerning courier delivery.) 5. For Wire Transfer Payment Procedures, please refer to the MDUFMA Fee Payment Instructions at the following URL: http://www.fda.gov/cdrh/mdufma/faqs.html#3a. You are responsible for paying all fees associated with wire transfers. 6. Include a copy of the completed Cover Sheet in volume one of the application when submitting to the FDA at either the CBER or CDRH Document Mail Center. 		
1. COMPANY NAME AND ADDRESS (Include name, street address, city, state, country, and post office code) SMITHS MEDICAL MD, INC. 1265 GREY FOX ROAD ARDEN HILLS, MN 55112 1.1 EMPLOYER IDENTIFICATION NUMBER (EIN) (b) (4)	2. CONTACT NAME DAVID SHORT 2.1 E-MAIL ADDRESS david.short@smiths-medical.com 2.2 TELEPHONE NUMBER (Include Area Code) 651-628-7240 2.3 FACSIMILE (FAX) NUMBER (Include Area Code) 651-628-7457	
3. TYPE OF PREMARKET APPLICATION (Select one of the following in each column; if you are unsure, please refer to the application descriptions at the following web site: http://www.fda.gov/oc/mdufma)		
Select an application type: <ul style="list-style-type: none"> <input checked="" type="checkbox"/> Premarket notification (510(k)); except for third party reviews <input type="checkbox"/> Biologics License Application (BLA) <input type="checkbox"/> Premarket Approval Application (PMA) <input type="checkbox"/> Modular PMA <input type="checkbox"/> Product Development Protocol (PDP) <input type="checkbox"/> Premarket Report (PMR) 		3.1 Select one of the types below: <ul style="list-style-type: none"> <input checked="" type="checkbox"/> Original Application Supplement Types: <ul style="list-style-type: none"> <input type="checkbox"/> Efficacy (BLA) <input type="checkbox"/> Panel Track (PMA, PMR, PDP) <input type="checkbox"/> Real-Time (PMA, PMR, PDP) <input type="checkbox"/> 180-day (PMA, PMR, PDP)
4. ARE YOU A SMALL BUSINESS? (See the instructions for more information on determining this status.)		
<input type="checkbox"/> YES, I meet the small business criteria and have submitted the required qualifying documents to FDA		<input checked="" type="checkbox"/> NO, I am not a small business
4.1 If Yes, please enter your Small Business Decision Number:		
5. IS THIS PREMARKET APPLICATION COVERED BY ANY OF THE FOLLOWING USER FEE EXCEPTIONS? IF SO, CHECK THE APPLICABLE EXCEPTION.		
<input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates, parents, and partner firms	<input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population	
<input type="checkbox"/> This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only	<input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially	
6. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION FOR WHICH FEES WERE WAIVED DUE TO SOLE USE IN A PEDIATRIC POPULATION THAT NOW PROPOSES CONDITION OF USE FOR ANY ADULT POPULATION? (If so, the application is subject to the fee that applies for an original premarket approval application (PMA).)		
<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO		
7. USER FEE PAYMENT AMOUNT SUBMITTED FOR THIS PREMARKET APPLICATION (FOR FISCAL YEAR 2004)		
(b) (4)		


Form FDA 3601 (08/2003)

DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION CDRH PREMARKET REVIEW SUBMISSION COVER SHEET		Form Approval OMB No. 9010-0120 Expiration Date: September 30, 2004. See OMB Statement on page 5	
Date of Submission 08/10/2004		User Fee Payment ID Number (b)(4) 3	
FDA Submission Document Number (if known)			
SECTION A TYPE OF SUBMISSION			
PMA <input type="checkbox"/> Original Submission <input type="checkbox"/> Premarket Report <input type="checkbox"/> Modular Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment <input type="checkbox"/> Licensing Agreement	PMA & HDE Supplement <input type="checkbox"/> Regular (120 day) <input type="checkbox"/> Special <input type="checkbox"/> Panel Track (PMA Only) <input type="checkbox"/> 30-day Supplement <input type="checkbox"/> 30-day Notice <input type="checkbox"/> 135-day Supplement <input type="checkbox"/> Real-time Review <input type="checkbox"/> Amendment to PMA & HDE Supplement <input type="checkbox"/> Other	PDP <input type="checkbox"/> Original PDP <input type="checkbox"/> Notice of Completion <input type="checkbox"/> Amendment to PDP	510(k) <input checked="" type="checkbox"/> Original Submission: <input type="checkbox"/> Traditional <input type="checkbox"/> Special <input type="checkbox"/> Abbreviated (Complete section I, Page 5) <input type="checkbox"/> Additional Information <input type="checkbox"/> Third Party
IDE <input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement		Humanitarian Device Exemption (HDE) <input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment	Class II Exemption Petition <input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information
Evaluation of Automatic Class III Designation (De Novo) <input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information		Meeting <input type="checkbox"/> Pre-510(K) Meeting <input type="checkbox"/> Pre-IDE Meeting <input type="checkbox"/> Pre-PMA Meeting <input type="checkbox"/> Pre-PDP Meeting <input type="checkbox"/> Day 100 Meeting <input type="checkbox"/> Agreement Meeting <input type="checkbox"/> Determination Meeting <input type="checkbox"/> Other (specify):	
Other Submission <input type="checkbox"/> 513(g) <input type="checkbox"/> Other (describe submission):			
Have you used or cited Standards in your submission? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No (If Yes, please complete Section I, Page 5)			
SECTION B SUBMITTER, APPLICANT OR SPONSOR			
Company / Institution Name Smiths Medical MD, Inc.		Establishment Registration Number (if known) 2183161	
Division Name (if applicable)		Phone Number (including area code) (651) 628-7240	
Street Address 1265 Grey Fox Road		FAX Number (including area code) (651) 628-4757	
City St. Paul	State / Province MN	ZIP/Postal Code 55112	Country USA
Contact Name David H. Short			
Contact Title Director Regulatory and Clinical Affairs		Contact E-mail Address david.short@smiths-medical.com	
SECTION C APPLICATION CORRESPONDENT (e.g., consultant, if different from above)			
Company / Institution Name			
Division Name (if applicable)		Phone Number (including area code) ()	
Street Address		FAX Number (including area code) ()	
City	State / Province	ZIP/Postal Code	Country
Contact Name			
Contact Title		Contact E-mail Address	

SECTION D1			REASON FOR APPLICATION - PMA, PDP, OR HDE		
<input type="checkbox"/> Withdrawal <input type="checkbox"/> Additional or Expanded Indications <input type="checkbox"/> Request for Extension <input type="checkbox"/> Post-approval Study Protocol <input type="checkbox"/> Request for Applicant Hold <input type="checkbox"/> Request for Removal of Applicant Hold <input type="checkbox"/> Request to Remove or Add Manufacturing Site	<input type="checkbox"/> Change in design, component, or specification: Software / Hardware Color Additive Material Specifications Other (specify below)	<input type="checkbox"/> Location change: Manufacturer Sterilizer Packager			
<input type="checkbox"/> Process change: Manufacturing Sterilization Packaging Other (specify below)	<input type="checkbox"/> Labeling change: Indications Instructions Performance Shelf Life Trade Name Other (specify below)	<input type="checkbox"/> Report Submission: Annual or Periodic Post-approval Study Adverse Reaction Device Defect Amendment			
<input type="checkbox"/> Response to FDA correspondence:		<input type="checkbox"/> Change in Ownership <input type="checkbox"/> Change in Correspondent <input type="checkbox"/> Change of Applicant Address			
<input type="checkbox"/> Other Reason (specify):					
SECTION D2			REASON FOR APPLICATION - IDE		
<input type="checkbox"/> New Device <input type="checkbox"/> New Indication <input type="checkbox"/> Addition of Institution <input type="checkbox"/> Expansion / Extension of Study <input type="checkbox"/> IRB Certification <input type="checkbox"/> Termination of Study <input type="checkbox"/> Withdrawal of Application <input type="checkbox"/> Unanticipated Adverse Effect <input type="checkbox"/> Notification of Emergency Use <input type="checkbox"/> Compassionate Use Request <input type="checkbox"/> Treatment IDE <input type="checkbox"/> Continued Access	<input type="checkbox"/> Change in: Correspondent / Applicant Design / Device Informed Consent Manufacturer Manufacturing Process Protocol - Feasibility Protocol - Other Sponsor	<input type="checkbox"/> Repose to FDA Letter Concerning: Conditional Approval Deemed Approved Deficient Final Report Deficient Progress Report Deficient investigator Report Disapproval Request Extension of Time to Respond to FDA Request Meeting Request Hearing			
<input type="checkbox"/> Report submission: Current Investigator Annual Progress Report Site Waiver Report Final					
<input type="checkbox"/> Other Reason (specify):					
SECTION D3			REASON FOR SUBMISSION - 510(k)		
<input checked="" type="checkbox"/> New Device	<input type="checkbox"/> Additional or Expanded indications	<input type="checkbox"/> Change in Technology			
<input type="checkbox"/> Other Reason (specify):					

SECTION E ADDITIONAL INFORMATION ON 510(K) SUBMISSIONS							
Product codes of devices to which substantial equivalence is claimed						Summary of, or statement concerning, safety and effectiveness information	
1	FPA	2		3		4	
5		6		7		8	
<input checked="" type="checkbox"/> 510 (k) summary attached <input type="checkbox"/> 510 (k) statement							
Information on devices to which substantial equivalence is claimed (if known)							
	510(k) Number		Trade or Proprietary or Model Name				Manufacturer
1	K974163	1	Sof-Set, model 315, 316	1			Medtronic MiniMed
2	K033892	2	UltraFlex Infusion Set	2			Disetronic Medical Systems
3	K972135	3	Comfort Set, formerly Maersk Medical Pureline Comfort SubQ Infusion Set	3			Unomedical A/S
4		4	NOTE: See TAB 13 for attached 510(k) Summaries.	4			
5		5		5			
6		6		6			
SECTION F PRODUCT INFORMATION - APPLICATION TO ALL APPLICATIONS							
Common or usual name or classification Intravascular Infusion Set							
	Trade or Proprietary or Model Name for This Device					Model Number	
1	Cleo 90 Infusion Set					1 21-7021, 21-7021, 21-7022, 21-7025, 21-7026,	
2						2 21-7027, 21-7030, 21-7031, 21-7032, 21-7035	
3						3 21-7036, 21-7037	
4						4	
5						5	
FDA document numbers of all prior related submissions (regardless of outcome)							
1	2	3	4	5	6	7	8
9	10	11	12				
Data included in Submission							
<input type="checkbox"/> Laboratory Testing <input type="checkbox"/> Animal Trials <input type="checkbox"/> Human Trials							
SECTION G PRODUCT CLASSIFICATION - APPLICATION TO ALL APPLICATIONS							
Product Code		C.F.R. Section (if applicable)			Device Class		
FPA		21 CFR 880.5440			<input type="checkbox"/> Class I <input checked="" type="checkbox"/> Class II <input type="checkbox"/> Class III <input type="checkbox"/> Unclassified		
Classification Panel 80 General Hospital							
Indications (from labeling) The Cleo 90 Infusion Set is designed for subcutaneous infusion of medication (including insulin) only. It consists of an applicator that contains an adhesive site with cannula and needle and a separate extension tubing set with buckle for connection to the site. After insertion of the needle and cannula, the applicator automatically retracts the needle, providing needle protection. The set is not intended for intravenous (IV) infusion or to infuse blood or blood products.							

<p><i>Note:</i> Submission of this information does not affect the need to submit a 2891 or 2891a Device Establishment Registration form.</p>		<p>FDA Document Number (if known)</p>	
SECTION H MANUFACTURING / PACKAGING / STERILIZATION SITES RELATING TO A SUBMISSION			
<input type="checkbox"/> Original <input checked="" type="checkbox"/> Add <input type="checkbox"/> Delete	FDA Establishment Registration Number (b)(4)	<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Manufacturer	<input checked="" type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager / Relabeler

<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete	FDA Establishment Registration Number	<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Manufacturer	<input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager / Relabeler
Company / Institution Name		Establishment Registration Number	
Division Name (if applicable)		Phone Number (including area code) ()	
Street Address		FAX Number (including area code) ()	
City	State / Province	ZIP Code	Country
Contact Name	Contact Title	Contact E-mail Address	

<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete	FDA Establishment Registration Number	<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Manufacturer	<input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager / Relabeler
Company / Institution Name		Establishment Registration Number	
Division Name (if applicable)		Phone Number (including area code) ()	
Street Address		FAX Number (including area code) ()	
City	State / Province	ZIP Code	Country
Contact Name	Contact Title	Contact E-mail Address	

SECTION I UTILIZATION OF STANDARDS					
Note: Complete this section if your application or submission cites standards or includes a "Declaration of Conformity to a Recognized Standard" statement.					
	Standards No.	Standards Organization	Standards Title	Version	Date
1	EN 1707	ISO:1997	Conical Fittings with a 6% (luer) taper for syringes, needles		
2	10555-5	ISO:1996	Sterile, single-use intravascular catheters- Part 5: Over-the-needle		
3	TIR 27:2001	AAMI	Sterilization of healthcare products - radiation sterilization		
4	11137	ASNI/AAMI:1994	Sterilization of Medical Devices ...		
5	11137-1	ASNI/AAMI:1994	Sterilization of Medical Devices ...		
6	11137-2	ASNI/AAMI:1994	Sterilization of Medical Devices ...		
7	11137-3	ASNI/AAMI	(Voting Draft) Sterilization of Medical Devices.....		
Please include any additional standards to be cited on a separate page.					
<p>Public reporting burden for this collection of information is estimated to average 0.5 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to:</p> <p style="text-align: center;">Food and Drug Administration CDRH (HFZ-342) 9200 Corporate Blvd. Rockville, MD 20850</p> <p><i>An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control</i></p>					

August 10, 2004

Smiths Medical MD, Inc.
Medication Delivery
and Patient Monitoring Division

1265 Grey Fox Road
St Paul MN 55112 USA
Tel: +1 651-633-2556 Fax: +1 651-628-7459
www.deltec.com

Food and Drug Administration
Office of Device Evaluation
Document Mail Center (HFZ-401)
Center for Devices and Radiological Health
9200 Corporate Blvd.
Rockville, MD 20850

RE: 510(k) Notification: Cleo™ 90 Infusion Set

Dear Sir or Madam:

This 510(k) Notification is to inform you of the intent of Smiths Medical MD, Inc. to market the substantially equivalent devices described below.

Provided in TAB 1 is a Table of Contents for this submission.

Provided in TAB 2 is a "Truthful and Accurate Statement."

Provided in TAB 3 is a 510(k) Summary of Safety and Effectiveness for the Cleo 90 Infusion Set.

Provided in TAB 4 is an Indications for Use statement for the Cleo 90 Infusion Set

Provided in TAB 5 is the Intravascular Administration Sets Checklist.

- 1) Classification name: Intravascular Administration Set
Common or usual name: Subcutaneous Infusion Set
Proprietary name: Cleo™ 90 Infusion Set
- 2) Establishment registration number: 2183502
- 3) Device classification: Class II per 21 CFR §880.5440, 80 FPA,
General Hospital and Personal Use Panel
- 4) Legally marketed device: N/A

- 5) Performance standards: Performance standards have not been promulgated to date for this device.
- 6) Proposed labeling: A sampling of proposed labeling for the Cleo 90 Infusion Set is enclosed within TAB 6. Also provided in TAB 6 is a listing of the device configurations that will be offered.
- No advertising has been developed as of the date of this submission.
- For comparative purposes, labeling for the Unomedical, Comfort Set is enclosed within TAB 7; the Disetronic, Ultraflex™ Set is enclosed within TAB 8; and the Medtronic MiniMed Sof-Set® is enclosed within TAB 9.
- 7) The Cleo 90 Infusion Set is similar to current commercially available subcutaneous infusion sets as follows:
- This set has the identical indication for use, see TAB 4.
 - The set is a similar design with a flexible cannula and site and a removable extension set with luer connection. See TAB 6 for configurations that will be offered.
 - The set is manufactured with similar materials. Provided in TAB 10 is a Biological Evaluation summary of the materials used in the Cleo 90 Infusion Set. All materials used and have patient or user contact have been found to be biocompatible and suitable for their intended use in the subject device.
 - Provided in TAB 11 is a summary of the packaging and sterilization information for the Cleo 90 Infusion Set.
- The Cleo 90 Infusion Set is different from current commercially available Cleo 90 Infusion Set as follows:
- The set incorporates a needle protection feature that is either activated automatically following insertion of the needle or by the user. See TAB 14 and 15 for additional information related to this feature.
 - Provided in TAB 12 is a description of the device and illustrations of the of the Cleo 90 Infusion Set.
- 8) Provided in TAB 13 is a comparison of the Cleo 90 Administration Set with other substantially equivalent devices that are currently commercially available. The chart provided highlights the similarities and differences of the devices compared.

- 9) Provided in TAB 14 is a summary of the *in-vitro* testing and a summary of the design control activities that details the verification and/or validation activities required. A Certification Statement stating compliance with these requirements is provided. Also provided in TAB 14 is the risk analysis conducted for the subject product.
- 10) The Cleo 90 Infusion Set needle incorporates an integral protection feature and is designed to reduce or eliminate needle sticks or injuries following use of the device. Evaluation of the needle protection feature was conducted in accordance with the "Supplementary Guidance on Premarket Notifications for Medical Devices with Sharps Injury Prevention Features; Guidance for industry and FDA", dated December 31, 2002.

Provided in TAB 15 is a copy of the protocol for the evaluation, including copies of the case report forms utilized; the final report of the findings from the evaluation; and copies of all completed evaluation forms from participants and observers who participated in the evaluation and that comprised the data used to evaluate the performance of the needle protection feature.

A total of (b) (4) individuals participated in the evaluation which included patients who routinely use an insulin pump and infusion set, parents or guardians of patients that insert an infusion set on behalf of a patient, and health care providers who instruct patients and caregivers on the use of infusion pumps and infusion sets. The evaluation demonstrated the Cleo 90 Infusion Set needle protection feature was safe and effective for its intended use. There were no needle sticks or injuries sustained by the participants during the evaluation.

If you have any questions or require additional information, please contact the undersigned.

Sincerely,
Smiths Medical MD, Inc.

David H. Short
Director Regulatory and Clinical Affairs
Tel. (651) 628-7240
Fax (651) 628-7457

510(k) Premarket Notification

Cleo™ 90 Infusion Set

Table of Contents

	TAB
Table of Contents	1
Truthful and Accurate Statement	2
510(k) Summary of Safety and Effectiveness for the Cleo™ 90 Infusion Set	3
Indications for Use Statement for the Cleo™ 90 Infusion Set	4
Intravascular Administration Sets Checklist	5
Proposed Labeling for the Cleo™ 90 Infusion Set	6
Cleo 90 Infusion Set Configurations	6
Cleo 90 Instructions for Use	6
Shelf Carton and Shelf Carton Label	6
End Cap Label for Applicator	6
Engraved Text for Applicator Housing	6
Unomedical Comfort Set Labeling	7
Disetronic Ultraflex™ Set Labeling	8
Medtronic MiniMed Soft-Set® Labeling	9
Biological Evaluation	10
Sterilization and Packaging Information	11
Summary of Sterilization Processes	11
Sterile Packaging Information	11
Final Report, Biological Challenge Test	11
Device Description and Illustrations	12
Cleo™ 90 Infusion Set, Exploded View	12
Dimensional Drawing Applicator	12
Substantially Equivalent Comparison	13
Comparison Chart of Cleo 90 Infusion Set	13
510(k) Summary Medtronic MiniMed Sof-Set®	13
510(k) Summary Unomedical Comfort Set	13
510(k) Summary Disetronic UltraFlex	13
<i>In-Vitro</i> Testing and Design Control	14
Summary of In-Vitro Testing	14
Design Control Activities	14
Certification Statement	14
Risk Management Summary	14
Simulated Clinical Use Evaluation	15

TAB 2

Truthful and Accurate Statement

**PREMARKET NOTIFICATION
TRUTHFUL AND ACCURATE STATEMENT
[As Required By 21 CFR 807.87(j)]**

I certify that, in my capacity as Director Regulatory and Clinical Affairs of Smiths Medical MD, Inc., I believe to the best of my knowledge, that all data and information submitted in the premarket notification are truthful and accurate and that no material fact has been omitted.

10 August 2007
Date

David H. Short
David H. Short
Director Regulatory and Clinical Affairs
Smiths Medical MD, Inc.

[Premarket Notification 510(k) Number]

[Cleo 90 Infusion Set]

TAB 3

**510(k) Summary of Safety and Effectiveness
for the Cleo™ 90 Infusion Set**

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

Cleo™ 90 Infusion Set

August 10, 2004

GENERAL INFORMATION

Applicant's Name and Address: Smiths Medical MD, Inc.
1265 Grey Fox Road
St. Paul, MN 55112

Contact Person: David H. Short
Director Regulatory and Clinical Affairs

Common/Usual Name: Subcutaneous Infusion Set

Proprietary Name: Cleo™ 90 Infusion Set

Equivalence Device Comparison: Unomedical Comfort Set; Disetronic® Ultraflex™ Set; and the Medtronic MiniMed Sof-Set®

II. DEVICE DESCRIPTION

The Cleo™ 90 Infusion Set consists of an applicator that contains an adhesive site with cannula and needle and separate extension tubing set with buckle. The infusion set has a standard luer connection. After insertion of the needle and cannula into the subcutaneous tissue, the applicator automatically retracts the needle into the applicator housing providing needle protection. The infusion set buckle attaches to the site allowing infusion of medication through the extension tubing.

III. INTENDED USE OF THE DEVICE

The Cleo 90 Infusion Set is designed for subcutaneous infusion of medication (including insulin) only. The set is not intended for intravenous (IV) infusion or to infuse blood or blood products.

IV. DEVICE COMPARISON

The Cleo 90 Infusion Set was compared to the following similar products: Unomedical, Comfort Set; Disetronic®, Ultraflex™ Set; and the Medtronic MiniMed, Sof-Set®. The Cleo 90 Infusion Set is substantially equivalent to these products by having the identical indication for use, a site with flexible catheter and needle for insertion into the subcutaneous tissue and a separate extension tubing with buckle for attachment to the site. The Unomedical Comfort Set utilizes a rigid needle rather than a flexible cannula. All sets have an adhesive patch that secures the catheter or needle to the skin following placement. The

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

Cleo™ 90 Infusion Set

August 10, 2004

Unomedical Comfort Set and Medtronic MiniMed Sof-Set have detachable extension sets, where as the Disetronic, Ultraflex Set has an integral extension set. For those devices with a flexible cannula, the needle is removed prior to infusion. The Cleo 90 Infusion Set incorporates a needle protection feature that retracts the needle automatically following insertion of the needle and cannula. The needle may also be retracted automatically by the user as well.

V. SUMMARY OF STUDIES

A. **Functional Testing**

In-vitro functional testing of the Cleo 90 Infusion Set was conducted. Biocompatibility testing was performed on the materials of the Cleo 90 Infusion Set.

B. **Clinical Studies**

Human clinical studies were not deemed necessary to evaluate the safety or effectiveness of the Cleo 90 Infusion Set.

To evaluate the safety and effectiveness of the needle protection feature of the Cleo 90 Infusion Set, a simulated clinical use evaluation was conducted in accordance with the "Supplementary Guidance on Premarket Notifications for Medical Devices with Sharps Injury Prevention Features; Guidance for Industry and FDA", December 31, 2002. Individuals who participated in the evaluation included patients who routinely use an insulin pump and infusion set, parents or guardians of patients who insert an infusion set on behalf of a patient, and health care providers who instruct patients and caregivers on the use of infusion pumps and infusion sets. This evaluation determined the needle protection feature operated as intended and no needle sticks or injuries were reported in the evaluation by the participants.

C. **Conclusions Drawn from the Studies**

The results of the testing conducted indicate the Cleo 90 Infusion Set functioned according to specifications and the materials used in the device are biocompatible. Based on these results, the product is considered acceptable for human use.

TAB 4

**Indications for Use Statement for the Cleo™ 90 Infusion
Set**

Indications for Use

510(k) Number (if known): _____

Device Name: Cleo 90 Infusion Set

Indications for Use:

The Cleo 90 Infusion Set is designed for subcutaneous infusion of medication (including insulin) only. The set is not intended for intravenous (IV) infusion or to infuse blood or blood products.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X OR Over-The Counter Use _____
(Per 21 CFR 801.109)

TAB 5

Intravascular Administration Sets Checklist

Intravascular Administration Set Checklist
Cleo™ 90 Infusion Set

Checklist Requirement	Section
1. Cover Letter	
2. Critical Elements	
Trade name	Cover Letter
Common name	Cover Letter
Classification name	Cover Letter
Procode/class	Cover Letter
Establishment registration number	Cover Letter
Purpose of submission	Cover Letter
Contact	Cover Letter
Telephone/fax	Cover Letter
3. Administration Requirements	
Indications for use	TAB 4 and 6
Truthful and accurate statement (21 CFR 807.87)	TAB 2
510(k) summary (21 CFR 807.92)	TAB 3
4. Labeling	
Identity of the device	Cover Letter, TAB 3 and 6
Intended use	TAB 3 and 6
Type, size/set length, ID, OD, drops/ml, priming volume, filter, needle gauge, length, quantity	TAB 6 and 13
Single use only	TAB 6
Sterile/non-toxic	TAB 6
Nonpyrogenic fluid path	TAB 6
Prescription statement	TAB 6
Frequency of replacement	TAB 6
Directions for use	TAB 6
Special limitations	TAB 6
Promotional materials	No advertising has been developed as of the date of this submission.
Plasticizer	TAB 10 (TOTM)
5. Description of Device	
Narrative description	Cover Letter, TAB 3 and 12
Photographs/drawing and/or labeled diagrams	TAB 12
Sample provided	Sample available upon request.
6. Intended Use	
Clear statement	TAB 4
7. Physical Specifications	
Dimensions and volumes	TAB 6 and 13
Proximal and distal end configuration	TAB 12 and 13
Connector type	TAB 12 and 13
Color (purpose)	Not applicable
Opacity	Not applicable
Markings and scales	Not applicable
Special features	Cover Letter and TAB 6
8. Mechanical Specifications	
Strength of materials	TAB 14
Stress characteristics	TAB 14
Fluid flow rate	TAB 14
Identify lubricant, if applicable	TAB 10
9. Material Specifications (Tests)	
Tensile strength	TAB 14
Burst test	TAB 14

Intravascular Administration Set Checklist
Cleo™ 90 Infusion Set

Checklist Requirement	Section
Leakage test	TAB 14
Security of attachments	TAB 14
Hardness	TAB 14
% Tubing elongation	TAB 14
Crimping/flexing with clamp	TAB 14
Cyclic performance with pump	Not applicable
10. Sterilization Information	
Method	TAB 11
Validation method	TAB 11
SAL	TAB 11
Packaging	TAB 11
EtO residuals	Not applicable
Pyrogen free method	TAB 11
Radiation dose	TAB 11
11. Biocompatibility Tests	
Cytotoxicity	TAB 10
Sensitization	TAB 10
Intracutaneous	TAB 10
Systemic toxicity	TAB 10
Hemocompatibility	TAB 10
Other	TAB 10
12. Descriptive Comparison to a Legally Marketed Device	
Identified appropriate legally marketed device(s)	Cover Letter, TAB 13
Side by side comparison	TAB 13
Features	Cover Letter, TAB 6 and 13
Intended use(s)	TAB 6 and 13
Labeling	TAB 6
All materials	TAB 10
Technological aspects	TAB 6, 12 and 13
Specifications	TAB 6 and 13
Performance	TAB 13
How differences may affect safety and effectiveness	Not applicable. Risk analysis conducted, TAB 14
13. Performance Data	
Risk analysis, if appropriate	TAB 14
Engineering	TAB 14
Microbial change, if appropriate	TAB 11
Simulated use study	TAB 14 and 15
Actual use study	Not applicable

TAB 6

Proposed Labeling for the Cleo™ 90 Infusion Set

Cleo 90 Infusion Set Configurations

Cleo 90 Instructions for Use

Shelf Carton and Shelf Carton Label

End Cap Label for Applicator

Engraved Text for Applicator Housing

Cleo 90 Infusion Set Configurations


Reorder Number (Variable Data)	Product Name	Catheter Length (Variable Data)	Tubing Length (Variable Data)	Number of Applicators (Variable Data)	Number of Extension Tubing Sets (Variable Data)
21-7220-24	Cleo 90	6 mm	24 inches (61 cm)	10	10
21-7221-24	Cleo 90	6 mm	31 inches (79 cm)	10	10
21-7222-24	Cleo 90	6 mm	42 inches (107 cm)	10	10
21-7225-24	Cleo 90	6 mm	24 inches (61 cm)	10	5
21-7226-24	Cleo 90	6 mm	31 inches (79 cm)	10	5
21-7227-24	Cleo 90	6 mm	42 inches (107 cm)	10	5
21-7230-24	Cleo 90	9 mm	24 inches (61 cm)	10	10
21-7231-24	Cleo 90	9 mm	31 inches (79 cm)	10	10
21-7232-24	Cleo 90	9 mm	42 inches (107 cm)	10	10
21-7235-24	Cleo 90	9 mm	24 inches (61 cm)	10	5
21-7236-24	Cleo 90	9 mm	31 inches (79 cm)	10	5
21-7237-24	Cleo 90	9 mm	42 inches (107 cm)	10	5

Note: Catheter Length, Tubing Length, Number of Applicators, and Number of Tubing Sets are variable data that is included in labeling when specific set is manufactured.

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

Fill the empty space in the cannula with the following prime:

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


ENGLISH

Indications for use:

The Comfort infusion set is indicated for subcutaneous infusion of medication, including insulin, administered by an insulin pump.

CONTRAINDICATIONS Comfort is indicated for subcutaneous use only and not for intravenous infusion or the infusion of blood or blood products.

WARNINGS

- The Comfort infusion set is only sterile and non-pyrogenic if the packaging is unopened and undamaged. Do not use it if the packaging is already open or is damaged.
- Go through instructions for use carefully before inserting Comfort, as failure to follow instructions may result in pain or injury.
- Inaccurate medication delivery, infection and/or site irritation may result from improper insertion and maintenance of infusion site.
- When using Comfort for the first time, do so in the presence of a healthcare provider.
- Comfort is a single use item, which must be disposed of after usage. Do not clean or re-sterilize.
- Be sure that the needle guard is removed before insertion.

- People with small or large amounts of subcutaneous fat should choose the insertion angle carefully, as the cannula could be placed in underlying muscle or dermal layer and hereby cause reduced or blocked medication delivery. Normal insertion angle between 20-45°. Consult your healthcare provider concerning this matter.
- Do not leave air in the infusion set. Make sure to fill Comfort completely.
- Do not in any way put disinfectants, perfumes, deodorants or other products containing alcohol or disinfectants in contact with the luer-lock connector and the tubing, as these may affect the integrity of the infusion set.
- Replace the infusion set if the adhesive tape becomes loose.
- Check the infusion set frequently to ensure that the soft cannula remains firmly in place, and replace it if it is not in place. Since the cannula is soft, it will not cause any pain if it the soft cannula becomes fully or partially dislodged from the skin, and this may therefore take place without notice. The soft cannula must always be completely inserted to receive full amount of medication.

- Replace infusion set every 48-72 hours according to Centers for Disease Control guidelines, or per healthcare provider's instructions.
- If infusion site becomes inflamed, replace set and use new site until the first site has healed.
- Do not re-insert introduced needle into infusion set. Reinsertion could cause tearing of soft cannula which would result in unpredictable medication flow.
- Never prime the infusion set or attempt to free a clogged tubing while set is inserted. You may accidentally inject too much medication.
- Use aseptic techniques when temporarily disconnecting Comfort, and seal the cannula housing with the cover provided. Consult with healthcare provider on how to compensate for missed medication when disconnected.
- Protect Comfort from direct sunlight. Store at room temperature.
- If infusing insulin, carefully monitor blood glucose levels when disconnected and after reconnecting.
- If infusing insulin, and blood glucose level becomes unexplainably high, or occlusion alarm occurs, check for clogs and leaks. If in doubt, change infusion set since soft cannula could be dislodged, crimped or partially clogged.

Discuss plan for rapid replacement of insulin with health care provider should any of these problems arise. Test blood glucose level to make sure problem is corrected.

- Check the blood glucose level 1-3 hours after introducing Comfort and check the blood glucose several times a day and measure your blood glucose regularly.
- If infusing insulin, do not change infusion set just prior to bedtime, unless blood glucose can be checked 1-3 hours after insertion.

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Deutsch

DISETRONIC

Manufactured for/
Manufacturer EU
Disetronic Medical Systems AG
Brunnmattstrasse 6
CH-3401 Burgdorf
Tel. +41 34 427 11 11
Fax +41 34 427 11 22
www.disetronic.com



Distributed by
Disetronic Medical Systems, Inc.
St. Paul, MN 55112-1014
Toll free 1-800-780-7801
www.disetronic-usa.com

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
Authorized Representative EU
Disetronic Medical Systems GmbH
DF-69366 Sulzbach/Ts

CE 0123


ID 420928/V03/08.01

Symbols on the label Símbolos en la etiqueta		Symboles de l'etiquette Symbole auf dem Etikett		REF	
	Consult instructions for use Leer las instrucciones para el uso Lire le mode d'emploi avant usage Gebrauchsanweisung beachten		Use by Fecha de caducidad Date limite d'utilisation Verwendbar bis		Reference number Número de artículo Numero d'article Bestellnummer
	Do not re-use No reutilizar Ne pas réutiliser Nicht zur Wiederverwendung		Date of manufacture Fecha de fabricación Date de fabrication Herstellungsdatum		Non Pyrogenic Apirógeno Apyrogène Pyrogenfrei
	Keep away from heat Proteger de calos Crainit la chaleur Vor Wärme schützen		Keep dry Proteger de humedad Crainit l'humidité Vor Feuchtigkeit schützen		PVC-free Sin PVC Sans PVC PVC-frei

INSTRUCTIONS FOR USE

Indications for use

Ultraflex™ is an infusion set for the subcutaneous infusion of insulin or medications administered with microdosage pumps (e. g. insulin pumps).

Contents of the Ultraflex™ II package

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WARNINGS

- This device is sterile unless the package is opened or damaged. Do not use if the sterile package has been opened or damaged.
- Protect the package from excessive humidity, sunlight and heat. Store at room temperature.
- Discard the set after each use. Do not clean or resterilize.
- The Ultraflex is for subcutaneous use only.
- Always tighten the luer lock firmly into the pump adapter. If this is not done, there is a danger of leakage of medication due to the excessively loose connection between the infusion set and the pump. Do not, however, use any tools to

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60

- Never use a soft cannula for more than 72 hours (3 days). Change the tubing when you change the cartridge; never use the tubing for more than 6 days.
 - If the infusion site becomes inflamed, replace the Ultraflex and use a new site.
 - Do not put disinfectants, perfumes, or deodorants on your Ultraflex as these may affect the integrity of the set.
 - Never prime the Ultraflex or attempt to free a clogged line while the infusion set is inserted. You may accidentally infuse an uncontrolled quantity of medication.
 - Avoid mechanical stresses on the site of infusion and on the Ultraflex. Pay attention when carrying heavy weights.
- Note for insulin use:**
- Check blood glucose level 1–3 hours after inserting the Ultraflex and monitor the infusion
- EN
- Check blood glucose frequently (check with your physician).
 - Do not insert a new soft cannula just prior to bedtime, unless blood glucose can be checked 1–3 hours after insertion.
 - If your blood glucose level becomes unexplainably high, or an «occlusion» alarm occurs, check for clogs and leaks; if in doubt, change Ultraflex: the soft cannula could be partially clogged. You must check at once whether the connection between the soft cannula and the cannula housing has separated. If this has become disconnected and has remained inside the body, you must seek medical assistance.
 - Check blood glucose 1–3 hours after changing Ultraflex and test again more often to make sure you have corrected the problem. If the blood glucose level is still high, contact your physician immediately. If your physician is not available, inject insulin manually as defined by your physician.
 - In case of nausea, bodily discomfort, vomiting, diarrhea, abdominal pain, strong thirst, apathy or tiredness, test your urine for ketones.

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Filling volume: U 100-Insulin	
12"/ 30 cm:	approx. 7 units approx. 70 µl
24"/ 60 cm:	approx. 10 units approx. 100 µl
31"/ 80 cm:	approx. 13 units approx. 130 µl
43"/ 110 cm:	approx. 18 units approx. 180 µl

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- Avoid hyperglycemia: You are personally responsible for this prevention. Together with your physician you should prepare a procedure for fast insulin replacement in the event that your soft cannula gets displaced. Interrupting the insulin supply can cause a life-threatening condition.
- Prepare the site to replace the soft cannula according to CDC (Centers for Disease Control and Prevention) guidelines or institutional protocol.

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9

Picture 7

To remove the introducer needle, hold the cannula housing firmly on your body. Place one finger on the top of the cannula housing (G) and withdraw the introducer guide needle in a straight manner using the needle handle (A).

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Important for users: It is advisable to administer a bolus of 0.2-0.5 unit of insulin U-100 (about 2-5 pl) in order to fill the empty space in the cannula after the introducer guide needle is withdrawn. Check with your health care team regarding the use of this procedure. Check blood glucose frequently.

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Picture 9

Hold adhesive (F) in place, gently press the side clips of the connector needle (H) and pull the connector needle (H) from the cannula housing (G).

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Biological Evaluation

TAB 11

Sterilization and Packaging Information

Summary of Sterilization Processes

Sterile Packaging Information

Final Report, Biological Challenge Test

101

TAB 12


Device Description and Illustrations

Cleo™ 90 Infusion Set, Exploded View




Dimensional Drawing Applicator


106

Cleo 90 Infusion Set Exploded View

107

Device State	Description of Device State
	<p>Cleo 90 Infusion Set applicator in shipped state, prior to removal of cap.</p>
	<p>Cleo 90 cap removed and applicator collar and spring fully compressed, exposing needle and cannula.</p>
	<p>Needle protection feature activated, showing needle fully retracted into housing. Site would be applied to skin, cannula inserted into subcutaneous tissue..</p>

Device State	Description of Device State
<p data-bbox="310 289 662 388">Needle on buckle pierces septum on site.</p>  <p data-bbox="472 573 646 646">Site</p>	<p data-bbox="954 310 1461 388">Extension tubing with buckle open and over site.</p> <p data-bbox="954 420 1388 493">Buckle can be placed in eight (8) separate orientations.</p> <p data-bbox="954 525 1494 598">Note: color of components is provided as an enhancement to the illustrations.</p>
<p data-bbox="332 777 673 871">Standard female luer fitting.</p> 	<p data-bbox="966 871 1469 945">Extension tubing with buckle showing buckle open with buckle on site.</p>
	<p data-bbox="974 1344 1485 1459">Extension tubing with buckle closed on site. Buckle needle will pierce site septum.</p>

TAB 13

Substantially Equivalent Comparison
Comparison Chart of Cleo 90 Infusion Set
510(k) Summary Medtronic MiniMed Sof-Set®
510(k) Summary Unomedical Comfort Set
510(k) Summary Disetronic UltraFlex

111

Comparison Table Substantial Equivalence

Set Name	Cleo 90	Sof-Set	Comfort	Ultraflex
Manufacturer	Smiths Medical	MiniMed	Unomedical	Disetronic
510(k) Number	-	K974163	K972135	K033892
Quick connect	At site	4" from site	At site	At site
Insertion	Manual	Auto/Manual	Manual	Manual
Needle protection	Auto/Cap	Silicone plug	Tube	Cap
Needle Removed from Cannula	Yes	Yes	No	Yes
Insertion Angle	90	90	20-45	90
Introducer Gauge	28	26	27	27
Cannula Lengths (mm)	6 & 9	6 & 9	17	8 & 10
Tubing Lengths (in)	24,31,42	24,42	24,31,43	24,31,43
Materials				
Tubing OD	PVC	PVC	PU	PU
Tubing ID	LDPE	HDPE	LDPE	PE
Luer	PC	PC	ABS	Poly (methacrylate)
Adhesive	UV	unknown	UV	unknown
Set male	PC	PVC	ABS	ABS
Set female	PC	PVC	na	na
Set needle	SST	SST	SST	SST
Site	PC	PC	PP	ABS
Site cover	PC	PC	Polyoxymethylene	na
Site septum	Silicone	Silicone	Silicone	Silicone
Skin Adhesive patch	Yes	No	Yes	Yes
Introducer needle	SST	SST	SST	SST
Introducer needle hub	PC	unknown	PP	unknown
Soft Cannula	FEP	FEP	PTFE	FEP
Insertion parts	PC	na	na	na
Insertion cap	HDPE	na	na	na

List of Abbreviations

Polyurethane	PU	Polyvinyl Chloride	PVC
Polyethylene	PE	Acrylonitrile, butadiene, styrene	ABS
High Density Polyethylene	HDPE	Fluorinated ethylene propylene	FEP
Low Density Polyethylene	LDPE	Polytetrafluoroethylene	PTFE
Polypropylene	PP	Stainless Steel	SST
Polycarbonate	PC	Not Applicable	na

FEB 2 1998

K974103

MiniMed Inc.
Premarket Notification - 510(k)
Sof-set Ultimate QR infusion set, models 315 and 316

Part E. 510(k) Summary

Submitter: MiniMed® Inc. 12744 San Fernando Rd., Sylmar, California 91342.

Contact: Don Selvey, Regulatory Affairs (818) 362-5958, 3011; (520) 527-0107 (v/f).

Name of Device: MiniMed® Sof-set® Ultimate QR® infusion sets, models 315 and 316.

Predicate Device: MiniMed Sof-set QR infusion sets, models 115 and 116.

Description of the New Device: The MiniMed Sof-set Ultimate QR infusion sets, models 315 and 316, are infusion administration sets, connecting to a medicine reservoir proximally, and inserted in the subcutaneous tissue of a user distally by means of an introducer needle. The reservoir to which the infusion set attaches proximally is inserted into an external infusion pump, such as the MiniMed 507 insulin pump.

The administration set attaches to the reservoir syringe by means of a female Luer connector, and subcutaneously in the user through an indwelling catheter made of FEP Teflon. The tubing is made of polyvinyl chloride (PVC) with a polyolefin liner. This configuration of PVC and polyolefin has been trademarked by MiniMed as Polyfin.

The 24 gauge indwelling catheter is introduced into the subcutaneous tissue by a removable 26 gauge introducer needle made of 304 stainless steel. The needle, indwelling catheter, and tubing share a common hub. The hub incorporates a winged configuration to facilitate handling of the administration set during insertion and stability following insertion.

On the skin-contact side of the wings, around the Teflon catheter, is an antibacterial dressing. Additionally, an adhesive dressing covers the wings of the administration set, securing the subcutaneous catheter and infusion line to the user.

Intended Use of the New Device: The MiniMed Sof-set Ultimate QR is intended for the subcutaneous infusion of medicine, including insulin, from an external infusion pump. The set is not intended nor indicated for use with blood or blood products.

©MiniMed, Sof-set, and QR are registered trademarks of MiniMed Inc.

000009

MiniMed Inc.
Premarket Notification - 510(k)
Sof-set Ultimate QR infusion set, models 315 and 316


Comparison of the Technological Features of the New Device and Predicate Device: The modified devices are substantially similar to the lawfully marketed predicate devices. Both are intended for subcutaneous delivery insulin or other appropriately labeled medicines from an infusion pump to the pump user. Both sets attach proximally to a reservoir by means of a Luer connector, and insert distally into the device user's subcutaneous tissues with a removable stainless steel needle and flexible catheter. Both devices are made of substantially similar materials.

Several differences exist between the new and predicate device:

- 1) changes in the shape/configuration of the hub and wings to be flatter and smoother;
- 2) the indwelling catheter is clear instead of white;
- 3) a change in the adhesive tab to be slightly longer and rotated 90°, facilitating tab removal prior to insertion and use with an auto inserter device;
- 4) a change in the adhesive tape: the infusion set will be secured with IV 3000, a product of Smith & Nephew;
- 5) a change in the polyethylene plastic in the inner layer of the tubing, from low density to a higher density. This has resulted in a more sensitive occlusion alarm, which will be promoted by MiniMed.

These modifications do not affect the safety or effectiveness of the device.

Signed,

Don Selvey
Regulatory Affairs
MiniMed Inc.

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FEB 2 1998

Mr. Don Selvey
Regulatory Affairs
MiniMed® Incorporated
12744 San Fernando Road
Sylmar, California 91342

Re: K974163
Trade Name: MiniMed® Sof-set® Ultimate QR® infusion
sets, models 315 and 316
Regulatory Class: II
Product Code: FPA
Dated: October 27, 1997
Received: November 5, 1997

Dear Mr. Selvey:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Pre-market Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does

115


Page 2 - Mr. Selvey

not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4618. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

Timothy A. Ulatowski
Director
Division of Dental, Infection Control
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

116

MiniMed Inc.
Premarket Notification - 510(k)
Sof-set Ultimate QR infusion set, models 315 and 316

INDICATIONS FOR USE

510(k) Number:

Device Name: MiniMed Sof-set Ultimate QR infusion sets, models 315 and 316.

Indications for Use: The MiniMed Sof-set Ultimate QR infusion set is indicated for the subcutaneous delivery of medicine, including insulin, from an external infusion pump.

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
Robert Curran
(Division Sign-Off)
Division of Dental, Infection Control,
and General Hospital Devices
510(k) Number 1979163

Prescription Use or Over-the-Counter Use
(Per 21 CFR 801.109)

000013

117

[FDA Home Page](#) | [CDRH Home Page](#) | [Search](#) | [CDRH A-Z Index](#) | [Contact CDRH](#)

[510\(k\)](#) | [Registration](#) | [Listing](#) | [Adverse Events](#) | [PMA](#) | [Classification](#) | [CLIA](#)
[CFR Title 21](#) | [Advisory Committees](#) | [Assembler](#) | [NHRIC](#) | [Guidance](#) | [Standards](#)

[New Search](#)

[Back To Search Results](#)

510(k) Premarket Notification Database

Device Classification Name	<u>Set, Administration, Intravascular</u>
510(K) Number	K974163
Regulation Number	<u>880.5440</u>
Device Name	MINIMED SOF-SET ULTIMATE QR INFUSION SET MODELS 3 <u>MEDTRONIC MINIMED</u>
Applicant	12744 San Fernando Rd. Sylmar, CA 91342
Contact	Don Selvey
Product Code	FPA
Date Received	11/05/1997
Decision Date	02/02/1998
Decision	Substantially Equivalent (SE)
Classification Advisory Committee	General Hospital
Review Advisory Committee	General Hospital
Statement/Summary/Purged Status	Summary/Purged 510(K)
Summary	<u>Summary</u>
Type	Traditional
Reviewed By Third Party	No
Expedited Review	No

Database Updated 8/05/2004

[CDRH Home Page](#) | [CDRH A-Z Index](#) | [Contact CDRH](#) | [Accessibility](#) | [Disclaimer](#)
[FDA Home Page](#) | [Search FDA Site](#) | [FDA A-Z Index](#) | [Contact FDA](#) | [HHS Home Page](#)

Center for Devices and Radiological Health / CDRH

file://F:\Projects\SQ_Infusion_Set\510(k)\510(k)_CoverSheet\k974163\K974163.htm

8/9/2004

118

510(K) SUMMARY

K972135

9.0 Summary of Safety and Effectiveness

AUG 26 1991

- 9.1 Trade/Proprietary Name: Maersk Medical Pureline™ Comfort™ Subcutaneous Infusion Set
- 9.2 Common/Usual Name: Subcutaneous Infusion Set with indwelling catheter
- 9.3 Classification Name: Intravascular Administration Set
- 9.4 Substantial Equivalence: The modified Pureline™ Comfort™ Subcutaneous Infusion Sets are substantially equivalent to the currently marketed Pureline™ Comfort™ Subcutaneous Infusion Sets (K95018).

9.5 Device Description

The general description, operation, construction and use of the Maersk Medical Pureline™ Comfort™ Subcutaneous Infusion Sets With Indwelling Catheters has not changed as a result of the modifications. They remain sterile, non-pyrogenic, single use subcutaneous infusion sets designed to be used with commercially available infusion devices. The following discussion remains valid.

There are two basic components provided for each device. The first is a stand-alone subcutaneous indwelling catheter. This catheter is provided as an integral assembly with a Teflon cannula, adhesive backed fixation tape, an injection port and the female portion of a proprietary plastic "click-lock" connector. The catheter assembly comes with a stainless steel insertion cannula. The insertion cannula is mounted to a male portion of the proprietary plastic "click-lock" connector. The insertion cannula is come to the user inserted through the injection port and the inner lumen of the Teflon catheter with the needle end protruding past the tip of the catheter. The male connector is locked to the female connector on the indwelling catheter. A needle protector is assembled over the Teflon catheter and the insertion cannula. A separate male portion of the proprietary connector without the insertion cannula is provided in the package. This component is used to attach to the female connector of the catheter after the indwelling Teflon cannula has been inserted and the steel insertion cannula has been withdrawn. The connector protects the indwelling catheter when the infusion set is not attached.

The second component of the Maersk Medical Subcutaneous infusion Set with Indwelling Catheter is the infusion set. The infusion set is comprised of a co-extruded tube with a female luer lock connector at the pump end and a stainless steel needle incorporated into the male portion of the proprietary plastic "click-lock" connector at the patient end. In order to maintain the sterility of the infusion path, a male luer cap covers the female luer lock connector and the male "click-lock" connector comes attached to a mating female connector. The sets come packaged in blister packs sealed with paper lid stock.

9.6 Intended Use

The Intended Use for the device has not changed. The Maersk Medical Pureline™ Comfort™ Subcutaneous Infusion Sets with Indwelling Catheters continue to be intended for the infusion and/or injection of fluids into the body below the surface of the skin. The indwelling catheter can be inserted

independently from the infusion catheter and can be accessed for injections through the injection port. The indwelling catheter can also be securely attached to the infusion catheter by means of a proprietary "click-lock" connector for the infusion of drugs subcutaneously. The infusion set can be detached from the indwelling catheter, and the catheter capped, to allow freedom from the infusion set and pump for showers, athletics or other activities.

9.7 Technological Characteristics

The technological characteristics of the device have not been affected by these modifications.

9.8 Performance Data

The results of the biocompatibility testing and component security testing showed the device is substantially equivalent to the unmodified device.

9.9 Conclusion

Based on the design equivalency and the functional and safety testing, Maersk Medical has determined that the Maersk Medical Pureline™ Comfort™ Indwelling Subcutaneous Catheter and Infusion Sets are substantially equivalent to the devices currently marketed in the United States.

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Mr. Lee Leichter
P/L Biomedical
7690 Cameron Circle
Fort Myers, Florida 33912

AUG 26 1997

Re: K972135
Trade Name: Maersk Medical Pureline Comfort Subcutaneous
Infusion Set
Regulatory Class: II
Product Code: FPA
Dated: May 14, 1997
Received: June 6, 1997

Dear Mr. Leichter:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Pre-market Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531


Page 2 - Mr. Leichter

through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4618. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

Timothy A. Ulatowski
Director
Division of Dental, Infection Control
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

INDICATIONS FOR USE

(510(k) File Number:

Device Name: Maersk Medical Pureline™ Comfort™ (Disetronic Tender™) Subcutaneous Infusion Set


Indications For Use: The Subcutaneous Infusion Sets with Indwelling Catheters are indicated for the infusion and/or injection of fluids into the body below the surface of the skin.

(Division Sign-Off) *Patricia Crockett*
Division of Dental, Infection Control,
and General Hospital Devices
510(k) Number K972135

Prescription Use
(Per 21 CFR 901.109)

PAGE 2

|23

[510\(k\)](#) | [Registration](#) | [Listing](#) | [Adverse Events](#) | [PMA](#) | [Classification](#) | [CLIA](#)
[CFR Title 21](#) | [Advisory Committees](#) | [Assembler](#) | [NHRIC](#) | [Guidance](#) | [Standards](#)

[New Search](#)

[Back To Search Results](#)

510(k) Premarket Notification Database

Device Classification Name	<u>Set, Administration, Intravascular</u>
510(K) Number	K972135
Regulation Number	880.5440
Device Name	MAERSK MEDICAL PURELINE COMFORT SUBCUTANEOUS INFU
Applicant	P/L BIOMEDICAL 7690 Cameron Circle Fort Myers, FL 33912
Product Code	FPA
Date Received	06/06/1997
Decision Date	08/26/1997
Decision	Substantially Equivalent (SE)
Classification Advisory Committee	General Hospital
Review Advisory Committee	General Hospital
Statement/Summary/Purged Status	Summary/Purged 510(K)
Summary	<u>Summary</u>
Type	Traditional
Reviewed By Third Party	No
Expedited Review	No

Database Updated 8/05/2004

[CDRH Home Page](#) | [CDRH A-Z Index](#) | [Contact CDRH](#) | [Accessibility](#) | [Disclaimer](#)
[FDA Home Page](#) | [Search FDA Site](#) | [FDA A-Z Index](#) | [Contact FDA](#) | [HHS Home Page](#)

Center for Devices and Radiological Health / CDRH

file://F:\Projects\SQ_Infusion_Set\510(k)\510(k)_CoverSheet\K972135\K972135.htm

8/9/2004

124

K033892

MAR - 9 2004

510(k) Summary

Introduction According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.

1) Submitter name, address, contact **Submitted by:**
 Disetronic Medical Systems AG
 Kirchbergstrasse 190, Postfach
 CH-3401 Burgdorf, Switzerland

United States Contact Person:
 Scott Thiel
 Roche Diagnostics
 9115 Hague Road
 Indianapolis, Indiana 46250
 317-521-3362
 scott.thiel@roche.com

Date Prepared: December 12, 2003

2) Device name Proprietary name: Ultraflex Infusion Set
 Common name: subcutaneous infusion set
 Classification name: intravascular administration set
 Product Code: FPA

3) Predicate device We claim substantial equivalence to the current legally marketed version of the same device.

4) Device Description The Ultraflex is a disconnectable infusion set with soft cannula perpendicular to the adhesive, for transfusion of insulin into the subcutaneous tissue. The unit is designed to interface with commercially available insulin infusion pumps with suitable connections. The insulin infusion pump systems are designed to control the delivery of insulin as prescribed by a health care professional. The system (infusion set, insulin infusion pump, and insulin) is indicated for patients with insulin dependent diabetes mellitus.

5) Intended use Ultraflex™ is an infusion set for the subcutaneous infusion of insulin administered with microdosage insulin pumps.

Continued on next page

510(k) Summary, Continued

6) Data demonstrating substantial equivalence

Testing of the modified Ultraflex Infusion Set demonstrated that the device meets the requirements for its intended use. The data also demonstrates that the Ultraflex Infusion Set is substantially equivalent to the predicate device.

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Re: K033892
Trade/Device Name: Ultraflex Infusion Set
Regulation Number: 880.5440
Regulation Name: Intravascular Administration Set
Regulatory Class: II
Product Code: FPA
Dated: December 12, 2003
Received: December 16, 2003

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If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4618. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,
Infection Control and Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

128

K033892

Indications for Use Statement

510(k) Number (if known):

Device Name: Ultraflex Infusion Set

Indications for Use:

Ultraflex™ is an infusion set for the subcutaneous infusion of insulin administered with microdosage insulin pumps.

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
[Back To Search Results](#)

510(k) Premarket Notification Database

Device Classification Name	<u>Set, Administration, Intravascular</u>
510(K) Number	<u>K033892</u>
Regulation Number	<u>880.5440</u>
Device Name	<u>ULTRAFLEX INFUSION SET</u> <u>DISETRONIC MEDICAL</u> <u>SYSTEMS AG</u>
Applicant	<u>Kirchberstrasse 190, Postfach</u> <u>Burgdorf, SW CH-34 01</u>
Contact	<u>Bernd Henningsen</u>
Product Code	<u>FPA</u>
Date Received	<u>12/16/2003</u>
Decision Date	<u>03/09/2004</u>
Decision	<u>Substantially Equivalent (SE)</u>
Classification Advisory Committee	<u>General Hospital</u>
Review Advisory Committee	<u>General Hospital</u>
Statement/Summary/Purged Status	<u>Summary Only</u>
Summary	<u>Summary</u>
Type	<u>Traditional</u>
Reviewed By Third Party	<u>No</u>
Expedited Review	<u>No</u>

Database Updated 8/05/2004

[CDRH Home Page](#) | [CDRH A-Z Index](#) | [Contact CDRH](#) | [Accessibility](#) | [Disclaimer](#)
[FDA Home Page](#) | [Search FDA Site](#) | [FDA A-Z Index](#) | [Contact FDA](#) | [HHS Home Page](#)

Center for Devices and Radiological Health / CDRH

file://F:\Projects\SQ_Infusion_Set\510(k)\510(k)_CoverSheet\k033892\K033892.htm

8/9/2004

130

TAB 14

In-Vitro Testing and Design Control

Summary of In-Vitro Testing

Design Control Activities

Certification Statement

Risk Management Summary

**Simulated Clinical Use Evaluation
of the
Smiths Medical MD, Inc.
Cleo™ 90 Infusion Set**

Final Report

August 6, 2004

**Smiths Medical MD, Inc.
1265 Grey Fox Road
St. Paul, MN 55112**

201

Infusion Set with subcutaneous catheter and tubing set

Use aseptic technique

Indications for Use - This infusions set is designed for *subcutaneous infusion of medication* (including insulin) *only*. It consists of an applicator that contains a self adhesive site with cannula and needle and a separate extension tubing set with buckle for connection to the site. After insertion of the needle and cannula, the applicator automatically retracts the needle, providing needle protection. The infusion set can be used with standard luer lock connectors.

Contraindications - **DO NOT** use for intravenous (IV) infusions, or to infuse blood or blood products.

Access Site Selection - Choose an area where you have a sufficient layer of subcutaneous tissue. Avoid areas that may be rubbed by clothing (waistband, etc.), and stay at least 2 inches (5cm) away from your belly button. Your healthcare provider can help you with site selection and give you information about properly caring for your site.

WARNINGS:

- Never prime (fill) the tubing with medication or insulin while the infusion set is connected to your body, or you could deliver an unintended bolus of medication.
- Ask your healthcare provider to identify proper technique and to help you select and maintain infusion sites:
 - Different people have different amounts of subcutaneous tissue. Ask your healthcare provider for help to identify the correct soft cannula length for you.
 - Improper insertion of the soft cannula and failure to maintain your insertion site and tubing can result in inaccurate medication delivery as well as irritation and/or infection of the infusion site.
 - Check your site regularly to make sure it is still in place. If the tape becomes loose, replace the infusion set as this may mean that the soft cannula is not completely inserted. If your infusion site becomes inflamed or infected, replace the infusion set and choose a new site well away from the previous site.
 - Remove all air in the tubing and cannula before starting infusion or medication delivery will be delayed as the tubing and cannula fill. Make sure the buckle is fully closed on the site hub or medication delivery will be interrupted.
 - Do not leave your tubing hanging loose where it could accidentally catch on objects and become disconnected.
- If using this infusion set for insulin delivery:
 - If your blood glucose becomes unexplainably high, or an occlusion alarm occurs, check the tubing and site - make sure the buckle is fully closed, check for clogs or leaks. If you cannot find the problem, replace your infusion set.
 - Develop a plan with your healthcare provider for quickly replacing missed insulin if you have a problem with the infusion set. You should also have a plan in place for when you temporarily disconnect your infusion set.
 - Test your blood glucose level regularly to quickly identify problems, and to make sure whatever treatment plan you are using has been effective. This is particularly important before you temporarily disconnect the tubing and when you reconnect.
 - Do not change your infusion set just before bedtime. You must be able to test your blood glucose 1 to 3 hours after insertion to make sure the infusion set is properly inserted and insulin delivery is working.

CAUTIONS:

- Applicator is **STERILE** and nonpyrogenic unless the tamper-evident ring is disconnected from the cap and/or it has been submerged in any liquid. Do not use if either of the above happens, or if applicator appears damaged or broken.
- Use applicator and tubing set once and then discard.
- Tubing set is **STERILE** and nonpyrogenic unless package is opened or damaged. Do not use if package is opened or damaged, or if product appears damaged.
- Do not put disinfectants, perfumes, deodorants or lotions on the infusion set or on the infusion site, as it could cause damage to the product.
- Store your infusion sets in a cool, dry place. Leaving infusion sets in direct sunlight or in your car may cause damage.
- Observe time and temperature limits for stability of medication. Make sure medication is at room temperature before use or air bubbles could form in tubing.
- Carefully read and follow instructions provided with pump, cartridge (if applicable) and medication.
- Dispose of used applicator, needles, tubing, cartridges (if applicable) and medication containers according to your local requirements.
- Do not use with substances that are incompatible with fluid path materials. Observe all time and temperature limits for drug/fluid stability.

Materials in the fluid path are medical grade: polycarbonate, silicone, fluorinated ethylene propylene, urethane/acrylate (UV cured), low-density polyethylene and stainless steel. Replace your infusion set every 48 to 72 hours, according to Centers for Disease Control guidelines, or as your healthcare provider instructs you to.

- continued -

CAUTION: Investigational device. Limited by Federal (United States) law to investigational use.

Symbols:

Latex-free

Expiration date

Gamma Sterilized

Reorder number

Single use


ATTENTION! See instructions for use

Lot No.

Manufactured in USA by or in Mexico for (See lot code for country of origin [3rd character from right]): A=USA, X=Mexico

2004-07

40-5374-01(03)

5374-01(03) IFU, Infusion Set

1

548

Disetronic® Ultraflex™

- English
- Español
- Français
- Deutsch

Instructions for use
Set for subcutaneous infusion

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Distributed by
Disetronic Medical Systems, Inc.
St. Paul, MN 55112-1014
Toll free 1-800-280-7801
www.disetronic-usa.com

Manufactured by Disetronic
in the USA

Authorized Representative EU
Disetronic Medical Systems GmbH
DE-95336 Sulzbach/Ts.

CE 0123

ID 420926/V03/08 01

	Symbols on the label Símbolos en la etiqueta	Symboles de l'étiquette Símbolos auf dem Etikett	REF
1		Consult instructions for use Leer las instrucciones para el uso Lire le mode d'emploi avant usage Gebrauchsanweisung beachten	
2		Use by Fecha de caducidad Date limite d'utilisation Verwendbar bis	
3		Do not re-use No reutilizar Ne pas réutiliser Nicht zur Wiederverwendung	
4		Keep away from heat Proteger de calor Craint la chaleur Vor Wärme schützen	
5			STERILE/O Sterilized using Ethylene Oxide Esterilizado con óxido de etileno Stérilisé à l'oxyde d'éthylène Sterilisiert mit Ethylenoxid
6			LOT Lot number Número de lote Número de lot Chargennummer
7			Non Pyrogenic Apirógeno Apyrogène Pyrogenfrei
8			PVC PVC-free Sin PVC Sans PVC PVC-frei

Indications for use
Ultraflex™ is an infusion set for the subcutaneous infusion of insulin or medications administered with microdosage pumps (e. g. insulin pumps).

Contents of the Ultraflex™ II package
inner packages 1 and 2.

Package 1 (complete set):

- A** introducer needle handle
- B** introducer needle protective cover
- C** introducer guide needle
- D** soft cannula
- E** adhesive paper backing
- F** self-adhesive
- G** cannula housing
- H** connector needle
- I** side clips of the connector needle
- J** tubing
- K** luer lock
- L** protective cap
- M** luer cap

Package 2 (additional base units):


Components A-G, I, from package 1

Contents of the Ultraflex™ I package

Package 1 only.

INSTRUCTIONS FOR USE

Instructions for use must be read carefully before use. Please also consult the manual of your insulin pump for information about pump therapy, sources of errors and potential risks and their causes.

Important: For first time use, set up must be performed in the presence of a health care professional.

WARNINGS

- This device is sterile unless the package is opened or damaged. Do not use if the sterile package has been opened or damaged.
- Protect the package from excessive humidity, sunlight and heat. Store at room temperature.
- Discard the set after each use. Do not clean or sterilize.
- The Ultraflex is for subcutaneous use only.
- Always tighten the luer lock firmly into the pump adapter. If this is not done, there is a danger of leakage of medication due to the excessively loose connection between the infusion set and the pump. Do not, however, use any tools to

screw the luer lock connector into place, since these could damage the material of the luer lock connector.

- Do not leave air bubbles in the Ultraflex. Prime completely.
- A small percentage of Ultraflex with soft cannulas may get crimped during insertion or displaced while in use. Should you note this problem please use a new Ultraflex as soon as possible.
- Check frequently to make sure the soft cannula remains firmly in place. Because the cannula is so soft, you will not notice any pain if it pulls out. The soft cannula must always be completely inserted, in order to receive the full amount of medication. Replace the Ultraflex if the dressing tape comes loose.
- Do not reinsert the introducer needle into the infusion set. Reinsertion could cause tearing of the soft cannula which would result in unpredictable medication flow.
- When temporarily disconnecting the Ultraflex pay attention to cleanliness of use. Always seal the tube with the enclosed protective cap.

- Never use a soft cannula for more than 72 hours (3 days). Change the tubing when you change the cartridge, never use the tubing for more than 6 days.
- If the infusion site becomes inflamed, replace the Ultraflex and use a new site.
- Do not put disinfectants, perfumes, or deodorants on your Ultraflex as these may affect the integrity of the set.
- Never prime the Ultraflex or attempt to free a clogged line while the infusion set is inserted. You may accidentally infuse an uncontrolled quantity of medication.
- Avoid mechanical stresses on the site of infusion and on the Ultraflex. Pay attention when carrying heavy weights.

Note for insulin use:

- Check blood glucose level 1–3 hours after inserting the Ultraflex and monitor the infusion

- site. Check blood glucose frequently (check with your physician).
- Do not insert a new soft cannula just prior to bedtime, unless blood glucose can be checked 1–3 hours after insertion.
- If your blood glucose level becomes unexplainably high, or an «occlusion» alarm occurs, check for clogs and leaks. If in doubt, change Ultraflex: the soft cannula could be partially clogged. You must check at once whether the connection between the soft cannula and the cannula housing has separated. If this has become disconnected and has remained inside the body, you must seek medical assistance.
- Check blood glucose 1–3 hours after changing Ultraflex and test again more often to make sure you have corrected the problem. If the blood glucose level is still high, contact your physician immediately. If your physician is not available, inject insulin manually as defined by your physician.
- In case of nausea, bodily discomfort, vomiting, diarrhea, abdominal pain, strong thirst, apathy or tiredness, test your urine for ketones.

5

6

- **Avoid hyperglycemia:** You are personally responsible for this prevention. Together with your physician you should prepare a procedure for soft cannula replacement in the event that your soft cannula gets displaced. Interrupting the insulin supply can cause a life-threatening condition.
- Prepare the site to replace the soft cannula according to CDC (Centers for Disease Control and Prevention) guidelines or institutional protocol.

manufacturer's user manual to prime the Ultraflex. Prime the tubing until medication drips from the needle tip. Do not leave air bubbles in either the cartridge or the tubing of the infusion set.

Filling volume: U 100-Insulin

12" / 30 cm:	approx. 7 units	approx. 70 µl
24" / 60 cm:	approx. 10 units	approx. 100 µl
31" / 80 cm:	approx. 13 units	approx. 130 µl
43" / 110 cm:	approx. 18 units	approx. 180 µl

Important: Do not use any tools to screw the luer lock connector into place. Tug on the tubing gently to ensure connection.

Change the cannula housing at least every 2-3 days (max. 72 hours) because of insulin resistance or risk of inflammation of the infusion site. Change the tubing when you change the cartridge. Never use the tubing for more than 6 days.

Insert Ultraflex™

Picture 1

Choose a site away from the belt line, bones, scar tissue, belly button and recent infusion sites. Clean hands and the site according to the instructions from your health care professional. Allow injection site to dry before the infusion set is inserted.

Picture 2

Remove luer lock cap (M). Attach the luer lock (K) to the pump by tightening it into the pump adapter. Follow the pump ma-

7

Picture 3

While grasping the cannula housing (G) of the Ultraflex by thumb and middle finger, remove the transparent protective cover (B) by pulling and turning. Keep the protective cover (B) for later use to dispose of the guide needle (C) in a safe and convenient manner.

Note: Keep the protective cap (L) in the individual blister pack. You will need this part to plug the tubing if you intend to dis-

8

connect from your pump.

Picture 4

Remove both pieces of the paper backing (E) from self-adhesive (F).

Picture 5

Pinch the subcutaneous tissue as you would to give an injection. Place your forefinger on top of the blue introducer needle (figure 5) and apply firm pressure when inserting the introducer needle (C) **quickly and straight** into the infusion site, at a 90° angle. Watch the cannula (D) while inserting to assure that it goes in smoothly without kinking or crimping.

Important: Make sure that the soft cannula is fully inserted into the subcutaneous tissue and not kinked at skin surface.

Note: Once the introducer needle punctures the skin, the guide needle and the soft cannula must be quickly pushed in the rest of the way, or the cannula may kink upon insertion, causing reduced or blocked medication delivery.

Picture 6

Carefully smooth out the bent self-adhesive (F) onto the skin and run your fingers over the self-adhesive dressing to ensure good skin contact is achieved.

9

Picture 7

To remove the introducer needle hold the cannula housing firmly on your body. Place one finger on the top of the cannula housing (G) and withdraw the introducer guide needle in a straight manner using the needle handle (A).

Note: The introducer guide needle of the Ultraflex goes through a white tube, called the soft cannula (D), which is what remains in the subcutaneous layer, after the introducer needle is removed.

Picture 8

Remember to dispose of the introducer needle in a safe manner (needle container). You may use the needle protective cover provided to cover the steel needle. The introducer guide needle and needle protective cover are made to lock together.

Important for users: It is advisable to administer a bolus of 0.2-0.5 unit of insulin U-100 (about 2.5µl) in order to fill the empty space in the cannula after the introducer guide needle is withdrawn. Check with your health care team regarding the use of this procedure. Check blood glucose frequently.

Disconnect Ultraflex™

Important for users: Please consult with your health care professional on how to compensate for any missed insulin during the time you are disconnected. Carefully

10

monitor blood glucose level when you are disconnected from the pump and after you reconnect.

Ultraflex allows you to temporarily disconnect from your pump for showers, swimming etc., without changing your Ultraflex.

Picture 9

Hold adhesive (F) in place, gently press the side clips of the connector needle (I) and pull the connector needle (H) from the cannula housing (G).

Picture 10

Insert protective cap (L) into the connector needle onto the recessed needle.

Picture 11

If desired, it is possible to position another protective cap on the cannula housing with the flat side facing upwards. A new protective cap is included in each single pack.

Reconnect Ultraflex™

Important: Before reconnecting, ensure that a drop of medication appears at the tip of the connector needle. Tug on the tubing gently to ensure connection after reconnecting. It is not necessary to give an additional bolus.

11

Picture 12

Remove the protective cap (L) and prime via pump. Fill the Ultraflex until medication comes out of the connecting needle. Start the pump afterwards.

If the protective cap was used on the cannula housing, hold adhesive (F) in place and remove protective cap.

Reattach the connector needle (H) to cannula housing (G). An audible «click» ensures the connector is locked.

Manufactured for

Disetronic Medical Systems AG
CH-3401 Burgdorf
Phone +41 34 427 11 11
Telefax +41 34 427 11 22
www.disetronic.com

Distributed by

Disetronic Medical Systems, Inc.
St. Paul, MN 55112-1014
Toll free 1-800-280-7801
www.disetronic-usa.com

Manufactured by Disetronic in the USA

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

Page 1 of 12

601

Photograph of Simulated Use "Skin" and Cleo™ Infusion Set and Disetronic® Ultraflex™

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Public Health Service
Food and Drug Administration
Memorandum

From: DATE: 20 Oct 2004
Reviewer(s) - Name(s) Scott A. Cousburn

Subject: 510(k) Number K042172

To: The Record - It is my recommendation that the subject 510(k) Notification:

- Refused to accept.
- Requires additional information (other than refuse to accept).
- Is substantially equivalent to marketed devices.
- NOT substantially equivalent to marketed devices.
- Other (e.g., exempt by regulation, not a device, duplicate, etc.)

- | | | |
|---|---|--|
| Is this device subject to Section 522 Postmarket Surveillance? | <input type="checkbox"/> YES | <input checked="" type="checkbox"/> NO |
| Is this device subject to the Tracking Regulation? | <input type="checkbox"/> YES | <input checked="" type="checkbox"/> NO |
| Was clinical data necessary to support the review of this 510(k)? | <input type="checkbox"/> YES | <input checked="" type="checkbox"/> NO |
| Is this a prescription device? | <input checked="" type="checkbox"/> YES | <input type="checkbox"/> NO |
| Was this 510(k) reviewed by a Third Party? | <input type="checkbox"/> YES | <input checked="" type="checkbox"/> NO |
| Special 510(k)? | <input type="checkbox"/> YES | <input checked="" type="checkbox"/> NO |
| Abbreviated 510(k)? Please fill out form on H Drive 510k/boilers | <input type="checkbox"/> YES | <input checked="" type="checkbox"/> NO |

- Truthful and Accurate Statement Requested Enclosed
- A 510(k) summary OR A 510(k) statement
- The required certification and summary for class III devices
- The indication for use form

Combination Product Category (Please see algorithm on H drive 510k/Boilers) N/A

Animal Tissue Source YES NO Material of Biological Origin YES NO

The submitter requests under 21 CFR 807.95 (doesn't apply for SEs):

- No Confidentiality
- Confidentiality for 90 days
- Continued Confidentiality exceeding 90 days

Predicate Product Code with class: Additional Product Code(s) with panel (optional):

Proc Code: FDA - Intravascular Administration SET
Class II 21 CFR 880.5440

Review: Michelle Hubbard for Anthony Watson H.O 10/7/04
(Branch Chief) (Branch Code) (Date)

Final Review: [Signature] 10/7/04
(Division Director) (Date)

REVISED: 3/14/95

THE 510(K) DOCUMENTATION FORMS ARE AVAILABLE ON THE LAN UNDER 510(K) BOILERPLATES TITLED "DOCUMENTATION" AND MUST BE FILLED OUT WITH EVERY FINAL DECISION (SE, NSE, NOT A DEVICE, ETC.).

"SUBSTANTIAL EQUIVALENCE" (SE) DECISION MAKING DOCUMENTATION

K042172

Reviewer: Scott Colburn

Division/Branch: DAGID/GHDB

Device Name: Product to Which Compared (510(K) Number If Known):

- Predicates: Sof-Set Infusion Set – MiniMed Inc. K974163;
- Maersk Medical Pureline Comfort Subcutaneous Infusion Set – Unomedical K0972135;
- Ultraflex Infusion Set - Disetronic Medical System AG K033892

YES NO

1. Is Product A Device	X		If NO = Stop
2. Is Device Subject To 510(k)?	X		If NO = Stop
3. Same Indication Statement?	X		If YES = Go To 5
4. Do Differences Alter The Effect Or Raise New Issues of Safety Or Effectiveness?			If YES = Stop NE
5. Same Technological Characteristics?	X		If YES = Go To 7
6. Could The New Characteristics Affect Safety Or Effectiveness?			If YES = Go To 8
7. Descriptive Characteristics Precise Enough?		X	If NO = Go To 10 If YES = Stop SE
8. New Types Of Safety Or Effectiveness Questions?			If YES = Stop NE
9. Accepted Scientific Methods Exist?			If NO = Stop NE
10. Performance Data Available?	X		If NO = Request Data
11. Data Demonstrate Equivalence?	X		Final Decision: SE

Note: In addition to completing the form on the LAN, "yes" responses to questions 4, 6, 8, and 11, and every "no" response requires an explanation.

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Device Description:

The Cleo™ 90 Infusion Set consists of an applicator that contains an adhesive site with cannula and needle and separate extension tubing set with buckle. The infusion set has a standard luer connection. After insertion of the needle and cannula into the subcutaneous tissue, the applicator automatically retracts the needle into the applicator housing providing needle protection. The infusion set buckle attached to the site allowing infusion of medication through the extension tubing.

Refer to SE 510 memo dated 07 Oct 2004

EXPLANATIONS TO "YES" AND "NO" ANSWERS TO QUESTIONS ON PAGE 1 AS NEEDED

7. Descriptive characteristics in the Cleo™ 90 Infusion Set alone would not be adequate to determine Substantial Equivalence. Simulated clinical use testing and performance testing to include; Pressure/Puncture tests, Tensile testing, and Mechanical testing was performed to determine safety and effectiveness. In addition, biocompatibility testing was done per ISO 10993.

11. Data presented in this file in the form of simulated clinical testing data, performance data, biocompatibility testing data demonstrate the Cleo™ 90 Infusion Set is substantially equivalent to its claimed predicates.

Page 1 of 510(k) review

**MEMO TO THE RECORD
510(K) REVIEW**

K042172

DATE: October 7, 2004
FROM: Scott Colburn BSN, RN

OFFICE: HFZ-480
DIVISION: DAGID/GHDB

COMPANY NAME: Smiths Medical MD, Inc.
DEVICE NAME: Cleo™ 90 Infusion Set

"SUBSTANTIAL EQUIVALENCE" (SE) DECISION-MAKING DOCUMENTATION

NARRATIVE DEVICE DESCRIPTION:

1. SUMMARY DESCRIPTION OF THE DEVICE UNDER REVIEW:
The Cleo™ 90 Infusion Set consists of an applicator that contains an adhesive site with cannula and needle and separate extension tubing set with buckle. The infusion set has a standard luer connection. The Infusion Set incorporates an integral protection feature and is designed to reduce needle sticks in the use of this device. After insertion of the needle and cannula into the subcutaneous tissue, the applicator automatically retracts the needle into the applicator housing providing needle protection. The infusion set buckle attached to the site allowing infusion of medication through the extension tubing. The buckle can be placed in eight separate orientations. The set is designed to be used with one hand for insertion of the device. However, individuals may vary in how they use the device and may use two hands without compromising safety or risk of needle stick (as demonstrated in simulated clinical study).
2. INTENDED USE:
The Cleo 90 Infusion Set is designed for subcutaneous infusion of medication (including insulin) only. The set is not intended for intravenous (IV) infusion or to infuse blood or blood products.
3. DEVICE DESCRIPTION:
 - A. Life-supporting or life-sustaining: No
 - B. Implant (short-term or long-term): No

8

Page 2 of 510(k) review

- C. Is the device sterile? Yes
If yes, is sterility information provided? Yes

(b) (4)

- D. Is the device for single use? Yes
- E. Is the device for prescription use? Yes
If yes, is prescription labeling included? Yes
- F. Is the device for home use or portable? Yes
- G. Does the device contain drug or biological product as a component? No
- H. Is this device a kit? Yes
If yes, and some or all of the components are not new, does the submission include a certification that these components were either preamendment or found to be substantially equivalent? No
- I. Software-driven: No
- J. Electrically Operated: No
- K. Applicable standards to which conformance has been demonstrated (e.g., IEC, ANSI, ASTM, etc.):
Standards to which testing was developed. Complete compliance of the following standards was not intended for this submission.
Applied Industrial Standards:
- ISO 10993-1 Biological evaluation of medical devices - Part 1: Evaluation and testing
 - AAMI TIR 21:2001
 - ISO 11137
 - ISO 11737-1, -2, -3
 - ISO 10555-1, -5
 - EN 1707
 - ISO 594-1, -2

Page 3 of 510(k) review

- L. Device(s) to which equivalence is claimed, manufacturer, and 510(k) number or preamendment status:
- Sof-Set Infusion Set – MiniMed Inc. K974163;
 - Maersk Medical Pureline Comfort Subcutaneous Infusion Set – Unomedical K0972135;
 - Ultraflex Infusion Set - Disetronic Medical System AG K033892

- M. Submission provides:
- | | |
|-------------------------------------|-----|
| a. comparative specifications | Yes |
| b. comparative in <u>vitro</u> data | No |
| c. performance data | Yes |
| d. animal testing | No |
| e. clinical testing | No |
| f. biocompatibility testing | Yes |

- N. Provide a statement of how the device is either similar to and/or different from other marketed devices, plus data (if necessary) to support the statement. Provide a summary about the devices design, materials, physical properties and toxicology profile if important.

The Cleo™ 90 Infusion Set was compared to three predicate devices; Unomedical comfort Set; Disetronic® Ultraflex™ Set; and the Medtronic MiniMed Sof-Set®. The Cleo™ 90 Infusion Set is similar in Indications for Use. The set is a similar design with a flexible cannula and adhesive site with a removable extension set with luer connection. Materials for main components are similar in nature and use and were testing for appropriate biocompatibility.

The Cleo™ 90 Infusion Set differs from its predicates by incorporating a needle protection/safety feature to reduce the risk of accidental needle sticks. This feature was tested under a simulated clinical study using participants who either were patients, parents or guardians of patients, and health care providers who instruct patients/guardians on the use of the infusion set. A total of 53 participants performing 530 insertion procedures were used in all. A majority of the participants (50) were able to use the device with one hand. However, some participants (3) felt more comfortable using both hands for placement. The study demonstrated the infusion set's safety feature was safe and effective under its intended use showing no needle sticks sustained or safety feature failures during the study.

Biocompatibility Testing:

Due to the multiple components which come into contact with the fluid path, testing was done in accordance with ISO 10993-1 for the following two areas; Surface device, skin, limited contact; and External Communicating device, Tissue/bone/dentin, permanent.

10

Page 4 of 510(k) review

Bench Testing:

Multiple bench performance tests were completed in the following areas of pressure/puncture, tensile, and mechanical tests as listed below:

(b)(4)

Based on the submission and additional information, I believe this device is substantially equivalent to the mentioned predicate devices stated in this submission. Sufficient information was submitted in a manner of simulated clinical, performance and engineering tests to substantiate the claim submitted for their device.


O. Does the submission include a summary of safety and effectiveness information upon which an equivalence determination is based? Yes

P. RECOMMENDATION:

I believe that this device is equivalent to: 80 FPA (Intravascular Administration Set)

Classification should be based on:

880.5440 Class: II

Scott Colburn BSN, RN

Neil Helbard
10/7/04

**SCREENING CHECKLIST
FOR ALL PREMARKET NOTIFICATION [510(k)] SUBMISSIONS**

510(k) Number: K042172

The cover letter clearly identifies the type of 510(k) submission as (Check the appropriate box):

- Special 510(k) - Do Sections 1 and 2
- Abbreviated 510(k) - Do Sections 1, 3 and 4
- Traditional 510(k) or no identification provided - Do Sections 1 and 4

Section 1: Required Elements for All Types of 510(k) submissions:

	Present or Adequate	Missing or Inadequate
Cover letter, containing the elements listed on page 3-2 of the Premarket Notification [510] Manual.	/	
Table of Contents.	/	
Truthful and Accurate Statement.	/	
Device's Trade Name, Device's Classification Name and Establishment Registration Number.	/	
Device Classification Regulation Number and Regulatory Status (Class I, Class II, Class III or Unclassified).	/	
Proposed Labeling including the material listed on page 3-4 of the Premarket Notification [510] Manual.	/	
Statement of Indications for Use that is on a separate page in the premarket submission.	/	
Substantial Equivalence Comparison, including comparisons of the new device with the predicate.	/	
510(k) Summary or 510(k) Statement.	/	
Description of the device (or modification of the device) including diagrams, engineering drawings, photographs or service manuals.	/	
Identification of legally marketed predicate device. *	/	
Compliance with performance standards. * [See Section 514 of the Act and 21 CFR 807.87 (d).]	/	
Class III Certification and Summary. **		
Financial Certification or Disclosure Statement for 510(k) notifications with a clinical study. * [See 21 CFR 807.87 (i)]		
510(k) Kit Certification ***		

- * - May not be applicable for Special 510(k)s.
- ** - Required for Class III devices, only.
- *** - See pages 3-12 and 3-13 in the Premarket Notification [510] Manual and the Convenience Kits Interim Regulatory Guidance.

Section 2: Required Elements for a SPECIAL 510(k) submission:

	Present	Inadequate or Missing
Name and 510(k) number of the submitter's own, unmodified predicate device.		
A description of the modified device and a comparison to the sponsor's predicate device.		
A statement that the intended use(s) and indications of the modified device, as described in its labeling are the same as the intended uses and indications for the submitter's unmodified predicate device.		
Reviewer's confirmation that the modification has not altered the fundamental scientific technology of the submitter's predicate device.		
A Design Control Activities Summary that includes the following elements (a-c):		
a. Identification of Risk Analysis method(s) used to assess the impact of the modification on the device and its components, and the results of the analysis.		
b. Based on the Risk Analysis, an identification of the required verification and validation activities, including the methods or tests used and the acceptance criteria to be applied.		
c. A Declaration of Conformity with design controls that includes the following statements:		
A statement that, as required by the risk analysis, all verification and validation activities were performed by the designated individual(s) and the results of the activities demonstrated that the predetermined acceptance criteria were met. This statement is signed by the individual responsible for those particular activities.		
A statement that the manufacturing facility is in conformance with the design control procedure requirements as specified in 21 CFR 820.30 and the records are available for review. This statement is signed by the individual responsible for those particular activities.		

Section 3: Required Elements for an ABBREVIATED 510(k)* submission:

	Present	Inadequate or Missing
For a submission, which relies on a guidance document and/or special control(s), a summary report that describes how the guidance and/or special control(s) was used to address the risks associated with the particular device type. (If a manufacturer elects to use an alternate approach to address a particular risk, sufficient detail should be provided to justify that approach.)		
For a submission, which relies on a recognized standard, a declaration of conformity [For a listing of the required elements of a declaration of conformity, SEE Required Elements for a Declaration of Conformity to a Recognized Standard, which is posted with the 510(k) boilers on the H drive.]		

For a submission, which relies on a recognized standard without a declaration of conformity, a statement that the manufacturer intends to conform to a recognized standard and that supporting data will be available before marketing the device.		
For a submission, which relies on a non-recognized standard that has been historically accepted by FDA, a statement that the manufacturer intends to conform to a recognized standard and that supporting data will be available before marketing the device.		
For a submission, which relies on a non-recognized standard that has <u>not</u> been historically accepted by FDA, a statement that the manufacturer intends to conform to a recognized standard and that supporting data will be available before marketing the device <u>and</u> any additional information requested by the reviewer in order to determine substantial equivalence.		
Any additional information, which is not covered by the guidance document, special control, recognized standard and/or non-recognized standard, in order to determine substantial equivalence.		

* - When completing the review of an abbreviated 510(k), please fill out an Abbreviated Standards Data Form (located on the H drive) and list all the guidance documents, special controls, recognized standards and/or non-recognized standards, which were noted by the sponsor.

Section 4: Additional Requirements for ABBREVIATED and TRADITIONAL 510(k) submissions (If Applicable):

	Present	Inadequate or Missing
a) Biocompatibility data for all patient-contacting materials, OR certification of identical material/formulation:	✓	
b) Sterilization and expiration dating information:	✓	
i) sterilization process	✓	
ii) validation method of sterilization process	✓	
iii) SAL	✓	
iv) packaging	✓	
v) specify pyrogen free	✓	
vi) ETO residues	✓	
vii) radiation dose	✓	
viii) Traditional Method or Non-Traditional Method		
c) Software Documentation:		

Items with checks in the "Present or Adequate" column do not require e additional information from the sponsor. Items with checks in the "Missing or Inadequate" column must be submitted before substantive review of the document.

Passed Screening Yes No

Reviewer: Scott A. Conner, D.D.


Concurrence by Review Branch: Wade Hubbard for Anthony Watson 10/7/04

Date: AUG 12 2004

Internal Administrative Form

	YES	NO
1. Did the firm request expedited review?		//
2. Did we grant expedited review?		//
3. Have you verified that the Document is labeled Class III for GMP purposes?		//
4. If, not, has POS been notified?		//
5. Is the product a device?	/	
6. Is the device exempt from 510(k) by regulation or policy?	/	/
7. Is the device subject to review by CDRH?	/	
8. Are you aware that this device has been the subject of a previous NSE decision?		
9. If yes, does this new 510(k) address the NSE issue(s), (e.g., performance data)?		//
10. Are you aware of the submitter being the subject of an integrity investigation?		/
11. If, yes, consult the ODE Integrity Officer.		
12. Has the ODE Integrity Officer given permission to proceed with the review? (Blue Book Memo #191-2 and Federal Register 90N0332, September 10, 1991).		/

510(k) "SUBSTANTIAL EQUIVALENCE" DECISION-MAKING PROCESS

- ❖ 510(k) Submissions compare new devices to marketed devices. FDA requests additional information if the relationship between marketed and "predicate" (pre-Amendments or reclassified post-Amendments) devices is unclear.
- ❖❖ This decision is normally based on descriptive information alone, but limited testing information is sometimes required.
- ❖❖❖ Data maybe in the 510(k), other 510(k)s, the Center's classification files, or the literature.

5

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- (GB) This product is covered by U.S. Patent 0,980,267 & 1,755,173, 1,176,662 & 2,323,380, 5,545,152 and 6,293,925. Other U.S. and/or foreign patents may be pending. MiniMed is a registered trademark of Medtronic MiniMed. Sof-set and Sof-serter are trademarks of Medtronic MiniMed.
(D) Sof-set and Sof-serter are covered by U.S. Patents 0,980,267 & 1,755,173, 1,176,662 & 2,323,380, 5,545,152 and 6,293,925. Other U.S. and/or foreign patents may be pending. MiniMed is a registered trademark of Medtronic MiniMed. Sof-set and Sof-serter are trademarks of Medtronic MiniMed.