

# U.S. Department of Health & Human Services

#### **Food and Drug Administration**

### **SAVE REQUEST**

USER: (lab)

**FOLDER:** K042841 - 365 pages

COMPANY: BIOMET, INC. (BIOMET)

PRODUCT: PROSTHESIS, HIP, SEMI-CONSTRAINED (METAL CEMENTED ACETABULAR COMPONENT) (JDL)

SUMMARY: Product: M2A/C2A ACETABULAR SYSTEM

DATE REQUESTED: Aug 11, 2014

**DATE PRINTED:** Aug 11, 2014

Note: Printed



DEC 21 2004



#### 510(k) Summary

Applicant/Sponsor:Biomet Manufacturing Corp.

Contact Person: Patricia Sandborn Beres

Senior Regulatory Specialist

Proprietary Name: M²a™/C²a™ Acetabular System

Common Name: Metallic Acetabular System

#### **Classification Name:**

1. Hip joint metal/metal semi-constrained, with an uncemented acetabular component prosthesis (21 CFR 888.3330)

2. Hip joint metal/metal semi-constrained, with a cemented acetabular component prosthesis (21 CFR 888.3320)

Legally Marketed Devices To Which Substantial Equivalence is Claimed: Biomet devices:

K993438 - Metal on Metal Acetabular System

K003363 - M²a™ 32mm Taper System

K861114 - Mallory/Head PF Acetabular Component

**Device Description**: The M²a™/C²a™ Acetabular System consists of a titanium outer acetabular shell with a cobalt alloy metallic liner for metal on metal articulation.

The acetabular shells are hemispherical in shape to closely match the natural acetabulum. Two screw holes in the dome allow for additional fixation by the use of screws. The outer surface of the shells are covered with Biomet's plasma sprayed coating.

The metallic cobalt alloy bearing liner fits into the outer shell by means of a taper similar to the taper used for the attachment of a modular head to a femoral stem. The metallic liners articulate with cobalt alloy modular heads.

MAILING ADDRESS P.O. Box 587 Warsaw, IN 46581-0587 SHIPPING ADDRESS 56 E. Bell Drive Warsaw, IN 46582

**OFFICE** 574.267.6639

FAX 574.267.8137 E-MAH. biomet@biomet.com

-K012841

M²a™/C²a™ Acetabular System 510(k) Summary Page 2

**Intended Use:** The M²a™/C²a™ Acetabular System is intended for cemented or non-cemented use in cases of:


- 1) Non-inflammatory degenerative joint disease including osteoarthritis and avascular necrosis
- 2) Rheumatoid arthritis
- 3) Correction of functional deformity
- 4) Revision procedures where other treatment or devices have been unsuccessful
- 5) Treatment of non-union, femoral neck fracture, trochanteric fractures of the proximal femur with head involvement, unmanageable using other techniques

**Summary of Technologies:** The technological characteristics of the new device are similar of identical to the predicates.

Non-Clinical Testing: None provided

Clinical Testing: None provided.

All trademarks are property of Biomet, Inc.





Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

DEC 21 2004

Ms. Patricia S. Andborn Beres Senior Regulatory Specialist 56 East Bell Drive P.O. Box 587 Warsaw, Indiana 46582

Re: K042841

Trade/Device Name: M<sup>2</sup>a / C<sup>2</sup>a<sup>TM</sup> Acetabular System

Regulation Number: 888.3330; 888.3320

Regulation Name: Hip joint metal/semi constrained, with uncemented acetabular component

prosthesis; Hip joint metal / metal semi-constrained acetabular

component prosthesis

Regulatory Class: III Product Code: KWA, JDL Dated: November 26, 2004 Received: November 29, 2004

Dear Ms. Andborn Beres:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0210. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <a href="http://www.fda.gov/cdrh/dsma/dsmamain.html">http://www.fda.gov/cdrh/dsma/dsmamain.html</a>

Sincerely yours,

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative and Neurological Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

## **Indications for Use**

510(k) Number (if known): とってる火	
Device Name: M²a™/C²a™ Acetabular Sys	stem_
Indications For Use:	
<ul> <li>cases of:</li> <li>1) Non-inflammatory degenerative joint dis necrosis</li> <li>2) Rheumatoid arthritis</li> <li>3) Correction of functional deformity</li> <li>4) Revision procedures where other treatments</li> </ul>	nent or devices have been unsuccessful acture, trochanteric fractures of the proximal
Prescription Use X AND/ (Part 21 CFR 801 Subpart D)	OR Over-The-Counter Use(21 CFR 807 Subpart C)
(PLEASE DO NOT WRITE BELOW THIS NEEDED)	LINE-CONTINUE ON ANOTHER PAGE IF
Concurrence of CDRH, Off	ice of Device Evaluation (ODE)
11.	Mar

Division of General, Restorative

and Neurological Devices

Page 1 of 1

510(k) Number K642841





Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

DEC 21 2004

Ms. Patricia S. Andborn Beres Senior Regulatory Specialist 56 East Bell Drive P.O. Box 587 Warsaw, Indiana 46582

Re: K042841

Trade/Device Name: M<sup>2</sup>a / C<sup>2</sup>a<sup>TM</sup> Acetabular System

Regulation Number: 888.3330; 888.3320

Regulation Name: Hip joint metal/semi constrained, with uncemented acetabular component

prosthesis; Hip joint metal / metal semi-constrained acetabular

component prosthesis

Regulatory Class: III Product Code: KWA, JDL Dated: November 26, 2004 Received: November 29, 2004

#### Dear Ms. Andborn Beres:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

#### Page 2 - Ms. Patricia S. Andborn Beres

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0210. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <a href="http://www.fda.gov/cdrh/dsma/dsmamain.html">http://www.fda.gov/cdrh/dsma/dsmamain.html</a>

Sincerely yours,

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative and Neurological Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

## **Indications for Use**

1110	210000000000000000000000000000000000000	
510(k) Number (if known): K of 2		
Device Name: M²a™/C²a™ Acetabula	ar Oyotom	
Indications For Use:		
<ul> <li>The M²a™/C²a™ Acetabular System is intended for cemented or non-cemented use in cases of:</li> <li>1) Non-inflammatory degenerative joint disease including osteoarthritis and avascular necrosis</li> <li>2) Rheumatoid arthritis</li> <li>3) Correction of functional deformity</li> <li>4) Revision procedures where other treatment or devices have been unsuccessful</li> </ul>		ling osteoarthritis and avascular ces have been unsuccessful
<ul><li>5) Treatment of non-union, femoral negative femur with head involvement, unm</li></ul>	ieck fracture, troc	Mantene fractures of the proximal
Prescription Use X (Part 21 CFR 801 Subpart D)	AND/OR	Over-The-Counter Use(21 CFR 807 Subpart C)
(PLEASE DO NOT WRITE BELOW NEEDED)	THIS LINE-CON	ITINUE ON ANOTHER PAGE IF

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)

Division of General, Restorative

and Neurological Devices

Page 1 of 1

510(k) Nouber K642841

Food and Drug Administration Center for Devices and Radiological Health Office of Device Evaluation Document Mail Center (HFZ-401) 9200 Corporate Blvd. Rockville, Maryland 20850

October 14, 2004

BIOMET, INC. 56 EAST BELL DR. P.O. BOX 587 WARSAW, IN 46582

ATTN: PATRICIA S. ANDBORN BERES

510(k) Number: K042841 Received: 14-0CT-2004

Product: M2A/C2A ACETABULAR

SYSTEM

The Food and Drug Administration (FDA), Center for Devices and Radiological Health (CDRH), has received the Premarket Notification you submitted in accordance with Section 510(k) of the Federal Food, Drug, and Cosmetic Act(Act) for the above referenced product. We have assigned your submission a unique 510(k) number that is cited above. Please refer prominently to this 510(k) number in any future correspondence that relates to this submission. We will notify you when the processing of your premarket notification has been completed or if any additional information is required. YOU MAY NOT PLACE THIS DEVICE INTO COMMERCIAL DISTRIBUTION UNTIL YOU RECEIVE A LETTER FROM FDA ALLOWING YOU TO DO SO.

On May 21, 2004, FDA issued a Guidance for Industry and FDA Staff entitled, "FDA and Industry Actions on Premarket Notification (510(k)) Submissions: Effect on FDA Review Clock and Performance Assessment". The purpose of this document is to assist agency staff and the device industry in understanding how various FDA and industry actions that may be taken on 510(k)s should affect the review clock for purposes of meeting the Medical Device User Fee and Modernization Act. Please review this document at http://www.fda.gov/cdrh/mdufma/guidance/1219.html

Please remember that all correspondence concerning your submission MUST be sent to the Document Mail Center (DMC)(HFZ-401) at the above letterhead address. Correspondence sent to any address other than the one above will not be considered as part of your official premarket notification submission. Also, please note the new Blue Book Memorandum regarding Fax and E-mail Policy entitled, "Fax and E-Mail Communication with Industry about Premarket Files Under Review". Please refer to this guidance for information on current fax and e-mail practices at www.fda.gov/cdrh/ode/a02-01.html.

You should be familiar with the regulatory requirements for medical device available at Device Advice http://www.fda.gov/cdrh/devadvice/". If you have other procedural or policy questions, or want information on how to check on the status of your submission, please contact DSMICA at (301) 443-6597 or its toll-free number (800) 638-2041, or at their Internet address http://www.fda.gov/cdrh/dsmamain.html or me at (301)594-1190.

Sincerely yours,

Marjorie Shulman Supervisory Consumer Safety Officer Office of Device Evaluation Center for Devices and Radiological Health

October 13, 2004

Document Mail Center (HFZ-401) Center for Devices and Radiological Health Food and Drug Administration 9200 Corporate Boulevard Rockville, Maryland 20850

RE:

M²a™/C²a™ Acetabular System

510(k) Premarket Notification Payment ID Number: 014934-956733

Dear Sir or Madam:

Enclosed is a 510(k) notification for the M²a™/C²a™ Acetabular System. We believe this device is substantially equivalent\* to other acetabular replacement devices on the market.

Please note, the C²a™ terminology in the title of this device refers to the C²a™-Taper Acetabular System that is the subject of a Modular Premarket Approval Application (PMA), M040011, submitted on September 29, 2004. The same acetabular shells are used for both the PMA ceramic-on-ceramic system and the metal-on-metal system which is the subject of this 510(k).

The sponsor of this 510(k) considers the existence of this notification confidential until a determination of substantial equivalence is made. Permission to fax or e-mail information related to this submission is granted by the Sponsor.

Sincerely,

Patricia Sandborn Beres

Senior Regulatory Specialist

aticia 5 Beres

Biomet Manufacturing Corp.

\*Any statement made in conjunction with this submission regarding and/or a determination of substantial equivalence to any other product is intended only to relate to whether the product can be lawfully marketed without pre-market approval or reclassification and is not intended to be interpreted as an admission or any other type of evidence in patent infringement litigation. [Establishment Registration and Premarket Notification Procedures. Final Regulation, Preamble, August 23, 1977, 42 FR 42520 (Docket No. 76N-0355)]

P.O. Box 587 Warsaw, IN 46581-0587

56 E. Bell Drive Warsaw, 1N 46582

OFFICE 574.267.6639

FAX 57/1.267.8137

E-MAIL biomet@biomet.com

#### **Table of Contents**

Medical Device User Fee Cover Sheet	İ
Premarket Notification Truthful and Accurate Statements2	2
Indications for Use Form	ŀ
<b>510(k) Notification</b> A. Administrative Information	5
B. Device Identification5	5
C. Device Descriptive Information6	)
D. 510(k) Summary	)
Exhibits Exhibit 1 Product Listing	2
Exhibit 2 Materials43	}
Exhibit 3 Sample Label 42 Package Insert 45 Surgical Technique 48	5
Exhibit 4 Predicate device comparison and information57	,
Exhibit 5 Testing125	;
Exhibit 6 510(k) Summary125	;
Appendix 1 Class III Certification and Summary244	ļ
Appendix 2 MAUDE Database Reports325	,

Form App	roved:OMB No. 0910-0511 Expiration Date: August 31, 2006. See instructions for OMB Statement.	
DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION PAYMENT	(b) (4)	
THE STATE OF THE PROPERTY OF T	TIDENTIFICATION NUMBI Payment Identification Num	
A completed Cover Sheet must accompany each original application or submit your application and fee payment:	applement subject to fees. The following actions must be taken to properly	
<ul> <li>Payment identification Number must be written on the check.</li> <li>Mail Check and Cover Sheet to the US Bank Lock Box, FDA Accepayment be submitted with the application.)</li> <li>If you prefer to send a check by a courier, the courier may deliver 956733, 1005 Convention Plaza, St. Louis, MO 63101. (Note: The if you have any questions concerning courier delivery.)</li> <li>For Wire Transfer Payment Procedures, please refer to the MDU http://www.fda.gov/cdrh/mdu/ma/fags.html#3a. You are responsible.</li> </ul>	k made payable to the Food and Drug Administration. Remember that the ount, P.O. Box 956733, St. Louis, MO 63195-6733. ( <i>Note: In no case should</i> the check and Cover Sheet to: US Bank, Attn: Government Lockbox is address is for courier delivery only. Contact the US Bank at 314-418-4821 FMA Fee Payment Instructions at the following URL:	
1. COMPANY NAME AND ADDRESS (Include name, street address, city state, country, and post office code)	, 2. CONTACT NAME PATRICIA BERES	
BIOMET MANUFACTURING CORP 56 EAST BELL DRIVE	2.1 E-MAIL ADDRESS patty.beres@biometmail.com	
P.O. BOX 587 WARSAW, IN 46581-0578	2.2 TELEPHONE NUMBER (Include Area Code) 574-267-6639	
1.1 EMPLOYER IDENTIFICATION NUMBER (EIN) 352074037	2.3 FACSIMILE (FAX) NUMBER (Include Area Code) 574-372-1683	
3. TYPE OF PREMARKET APPLICATION (Select one of the following in at the following web site: http://www.fda.gov/oc/mdufma	each column; if you are unsure, please refer to the application descriptions	
Select an application type:	3.1 Select one of the types below:	
Premarket notification (510(k)); except for third party reviews	☑ Original Application	
Biologics License Application (BLA)	Supplement Types:	
Premarket Approval Application (PMA)	Efficacy (BLA)	
☐ Modular PMA	Panel Track (PMA, PMR, PDP)	
Product Development Protocol (PDP)	Real-Time (PMA, PMR, PDP)	
Premarket Report (PMR)		
4. ARE YOU A SMALL BUSINESS? (See the instructions for more inform	ation on determining this status.)	
YES, I meet the small business criteria and have submitted the required qualifying documents to FDA	NO, I am not a small business	
4.1 If Yes, please enter your Small Business Decision Number:		
5. IS THIS PREMARKET APPLICATION COVERED BY ANY OF THE FOREXCEPTION.	DLLOWING USER FEE EXCEPTIONS? IF SO, CHECK THE APPLICABLE	
This application is the first PMA submitted by a qualified small busin including any affiliates, parents, and partner firms	ness,	
This biologics application is submitted under section 351 of the Pub Health Service Act for a product licensed for further manufacturing only		
6. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION FOR W POPULATION THAT NOW PROPOSES CONDITION OF USE FOR ANY applies for an original premarket approval application (PMA).)		
☐ YES ☑ NO		
7. USER FEE PAYMENT AMOUNT SUBMITTED FOR THIS PREMARK	ET APPLICATION (FOR FISCAL YEAR 2004)	

## PREMARKET NOTIFICATION TRUTHFUL AND ACCURATE STATEMENT

(As Required by 21 CFR 807.87(j))

I certify, in my capacity as a Development Engineer of Biomet Manufacturing Corp., I believe to the best of my knowledge, that all data and information submitted in the premarket notification are truthful and accurate and that no material fact has been omitted.

Signature

W. Jason Slone

Typed Name

Date

M²a™/C²a™ Acetabular System

Device

## PREMARKET NOTIFICATION TRUTHFUL AND ACCURATE STATEMENT

(As Required by 21 CFR 807.87(j))

I certify, in my capacity as Vice President of Regulatory Affairs and Quality Assurance, Biomet Inc., I believe to the best of my knowledge, that all data and information submitted in the premarket notification are truthful and accurate and that no material fact has been omitted.

Robert E. Durgin

Typed Name

/2 OcT 0 4 Date

M²a™/C²a™ Acetabular System

Device

## **Indications for Use**

510(k) Number (if known):		
Device Name: M²a™/C²a™ Acetab	ular System	
Indications For Use:		
<ul> <li>The M²a™/C²a™_Acetabular System cases of:</li> <li>1) Non-inflammatory degenerative junctosis</li> <li>2) Rheumatoid arthritis</li> <li>3) Correction of functional deformity</li> <li>4) Revision procedures where othe</li> <li>5) Treatment of non-union, femoral femur with head involvement, un</li> </ul>	oint disease inclu / r treatment or de neck fracture, tro	uding osteoarthritis and avascular vices have been unsuccessful ochanteric fractures of the proximal
Prescription Use X (Part 21 CFR 801 Subpart D)	AND/OR	Over-The-Counter Use (21 CFR 807 Subpart C)
(PLEASE DO NOT WRITE BELOW NEEDED)	THIS LINE-COI	NTINUE ON ANOTHER PAGE IF
Concurrence of CDF	RH, Office of Dev	rice Evaluation (ODE)

Page 1 of 1

#### 510(k) Notification

#### A. ADMINISTRATIVE INFORMATION

**Applicant or Sponsor:** Biomet Manufacturing Corp.

56 East Bell Drive P.O. Box 587

Warsaw, Indiana 46581-0587

Contact Person: Patricia Sandborn Beres

Senior Regulatory Specialist Biomet Manufacturing Corp.

P.O. Box 587

Warsaw, Indiana 46581-0587

Phone: (574) 267-6639 FAX: (574) 372-1683

E-Mail: patty.beres@biometmail.com

#### Manufacturing Site(s):

Specification holder:

Biomet Manufacturing Corp.

56 East Bell Drive

Warsaw, Indiana 46582

Establishment Registration Number: 1825034

Manufacturer/Contract Manufacturer:

Biomet Manufacturing Corp.

56 East Bell Drive

Warsaw, Indiana 46582

Establishment Registration Number: 1825034

(4) Contract Storilizor(s):

#### B. DEVICE IDENTIFICATION

Proprietary Name: M²a™/C²a™ Acetabular System

Common or Usual Name: Metallic Acetabular System

#### Classification Name:

- 1. Hip joint metal/metal semi-constrained, with an uncemented acetabular component prosthesis (21 CFR 888.3330)
- 2. Hip joint metal/metal semi-constrained, with a cemented acetabular component prosthesis (21 CFR 888.3320)

#### **Device Classification:**

- 1. Preamendment Class III
- 2. Preadmendment Class III

A Class III Certification and Summary may be found in Appendix 1 including the Medical Device Reports tabular summary and copies of the MAUDE Database reports. Class III literature review and copies of the abstracts are located in Appendix 2.

#### **Device Product Code:**

- 1. 87 KWA
- 2. 87 JDL

**Performance Standards/Guidance Documents:** No performance standards have been developed for this type of device.

Previous FDA Status: Components of this system were previously cleared through 510(k) K993438 - Metal on Metal Acetabular System and K003363 - M²a™ 32mm Taper System for non-cemented application only. The shell of this system is currently the subject of a PMA application for ceramic-on-ceramic articulation (M040011).

#### C. DEVICE DESCRIPTIVE INFORMATION

**Intended Use:** The M²a™/C²a™ Acetabular System is intended for cemented or uncemented use in cases of:

- 1) Non-inflammatory degenerative joint disease including osteoarthritis and avascular necrosis
- 2) Rheumatoid arthritis
- 3) Correction of functional deformity
- 4) Revision procedures where other treatment or devices have been unsuccessful
- 5) Treatment of non-union, femoral neck fracture, trochanteric fractures of the proximal femur with head involvement, unmanageable using other techniques

**Device Description:** The M²a™/C²a™ Acetabular System consists of a titanium outer acetabular shell with a cobalt alloy metallic liner for metal on metal articulation. A cobalt alloy modular head completes the system.

The acetabular shells are hemispherical in shape to closely match the natural acetabulum. The shells are available in outer diameters of 48mm to 70mm in 2 mm increments. The shell features eight radial fins to aid in the prevention of rotation. Two screw holes in the dome allow for additional fixation by the use of 6.5mm screws. The outer surface of the shells are covered with Biomet's plasma spray coating.

The metallic liners contained in this submission are identical to those previously cleared in 510(k) K993438 - Metal on Metal Acetabular System and K003363 - M²a™ 32mm Taper System. The taper locking mechanism of the shells is also identical to that of the shells cleared in these two 510(k)s.

The metallic, cobalt alloy bearing liner fits into the outer shell by means of a taper similar to the taper used for the attachment of a modular head to a femoral stem.

The locking mechanism consists of an 18°55' taper (b) (4)

(b) (4)

(b) (4)

I his is less than half the allowable deviation set forth in the ISO standard for modular head sphericity (ISO 7206-2).

The metallic liners articulate with cobalt alloy modular heads identical to those cleared through 510(k) K993438 - Metal on Metal Acetabular System and K003363 - M²a™ 32mm Taper System. (b) (4)

modular heads are available in seven neck lengths ranging from -6mm to +12mm.

Each modular head has Biomet's Type I Taper and will mate with any Biomet Type I Taper femoral component.

A part number listing and engineering drawings can be found in **Exhibit 1**. In addition, there are manual surgical instruments listed that are used with the device.

Materials: The acetabular shell is manufactured from titanium alloy (Ti-6Al-4V) conforming to ASTM F-136 or F-1472. The outer surface of the shell is plasma spray coated with Ti-6Al-4V powder conforming to ASTM F-1580. Please see the sample engineering drawings contained in Exhibit 1 for the location of the coating. A full characterization of this coating may be found in Master File MAF-153. A table summarizing this information is contained in Exhibit 2. This coating is identical to the coatings cleared in 510(k) K993438 - Metal on Metal Acetabular System and K003363 - M²a™ 32mm Taper System.

The metallic acetabular liner and metallic modular heads are manufactured from cobalt alloy (Co-Cr-Mo) conforming to ASTM F-1537.

**Labeling:** Copies of the package label, package insert, surgical technique may be found in **Exhibit 3**.

**Sterility Information:** Devices are provided sterile by radiation methods as follows:

Radiation Type: Gamma
Radiation Source: Cobalt 60
Minimum Dosage: 25 kGy
Maximum Dosage: 40 kGy
Sterility Assurance Level: 10-6

Sterility Validation Method: AMMI/ISO 11137, Method 1

Pyrogen-Free: no claims will be made

Labeling: All packages will display a yellow to red chemical indication dot along with a statement that the device has been sterilized by gamma irradiation, 2.5 Mrads.

Sterilization Site:



A summary of Biomet's sterilization methods is presented in Biomet's Masterfile MAF-153. Sample bioburden audits and pyrogen testing is also included. Both of these tests are done on a periodic basis for all Biomet devices.

**Packaging Description**: Each component is supplied in an individual sterile package. They are then placed in an inner blister pack sealed with a Tyvek® lid which fits into an outer blister pack also sealed with a Tyvek® lid. The entire unit is placed in a cardboard box, shrink wrapped for protection.

Substantial Equivalence: Exhibit 4 contains a table comparing the M²a™/C²a™ Acetabular System to several of Biomet's cleared devices. These devices are:

K993438 - Metal on Metal Acetabular System

K003363 - M2a 32mm Taper System

K861114 - Mallory/Head PF Acetabular Component

#### Indications for Use

The M²a™/C²a™ Acetabular System has indications for use similar to the predicates. Predicate Mallory/Head acetabular shells have been cleared for both cemented and non-cemented applications.

#### **Technological Characteristics**

The acetabular shell diameters and overall geometry are identical to the Mallory/Head shells cleared in 510(k)s K993438 and K003363. The only difference is the addition of 2 screw holes to the Mallory/Head configuration. Besides the Mallory/Head style shells, the predicate 510(k)s contained the Universal style acetabular shell with 2 screw holes in the dome.

**8** 2 *&* 

The shells inner taper angle that provides fixation of the metallic liner is identical to the predicates. Slight dimensional changes have been made to the inner shell geometry to insure proper seating of the liner in the shell even when the tolerances are at their limits. The surface roughness of the shell's inner taper has been increased over the predicate. In order to confirm that the locking strength of the components has not been altered from the predicates, testing has been conducted. **Exhibit 5** contains the test report and a discussion to the relevance of this testing to the current device.

The metallic liners are identical to those cleared in 510(k) s K993438 and K003363.

#### Summary

In summary, the M²a™/C²a™ Acetabular System is substantially equivalent to the predicate device in that:

- They have similar indications for use
- The acetabular shells are identical in overall design with no new features
- The taper locking mechanism for the metal liner is identical

#### D. 510(k) SUMMARY

A summary of information pertaining to the safety and effectiveness of this type of device is contained in Exhibit 6.

9 29

All trademarks are owned by Biomet, Inc .except for the following: Tyvek is a trademark of trademark of E.I. duPont de Nemours and Company

### **DEVICE LISTING**

## **Liner Component Listing**

Part Number	Description
15-105000 15-105002	28mm I.D./37 Taper 28mm I.D./41 Taper
15-105004	32mm I.D./41 Taper

## 28mm Modular Head Component Listing

Part Number	Description
11-163660	Metal on Metal 28mm -6mm Modular Head
11-163661	Metal on Metal 28mm -3mm Modular Head
11-163662	Metal on Metal 28mm std. Modular Head
11-163663	Metal on Metal 28mm +3mm Modular Head
11-163664	Metal on Metal 28mm +6mm Modular Head
11-163665	Metal on Metal 28mm +9mm Modular Head
11-163666	Metal on Metal 28mm +12mm Modular Head

## 32mm Modular Head Component Listing

Part Number	Description
11-163667	M²a™ 32mm Modular Head –6mm
11-163668	M²a™ 32mm Modular Head –3mm
11-163669	M²a™ 32mm Modular Head std
11-163670	M²a™ 32mm Modular Head +3mm
11-163671	M²a™ 32mm Modular Head +6mm
11-163672	M²a™ 32mm Modular Head +9mm skirt
11-163673	M²a™ 32mm Modular Head +12mm skirt

## **Acetabular Shell Component Listing**

Part Number	Description
10-111148	37mm x 48mm Mallory-Head® Radial 2-Hole Shells
10-111150	37mm x 50mm Mallory-Head® Radial 2-Hole Shells
10-111152	41mm x 52mm Mallory-Head® Radial 2-Hole Shells
10-111154	41mm x 54mm Mallory-Head® Radial 2-Hole Shells
10-111156	41mm x 56mm Mallory-Head® Radial 2-Hole Shells
10-111158	41mm x 58mm Mallory-Head® Radial 2-Hole Shells
10-111160	41mm x 60mm Mallory-Head® Radial 2-Hole Shells
10-111162	41mm x 62mm Mallory-Head® Radial 2-Hole Shells
10-111164	41mm x 64mm Mallory-Head® Radial 2-Hole Shells
10-111166	41mm x 66mm Mallory-Head® Radial 2-Hole Shells
10-111168	41mm x 68mm Mallory-Head® Radial 2-Hole Shells
10-111170	41mm x 70mm Mallory-Head® Radial 2-Hole Shells


## **Screw Component Listing**


Part Number	Description
103530	Titanium Acetabular Screw 6.5 x 15mm
103531	Titanium Acetabular Screw 6.5 x 20mm
103532	Titanium Acetabular Screw 6.5 x 25mm
103533	Titanium Acetabular Screw 6.5 x 30mm
103534	Titanium Acetabular Screw 6.5 x 35mm
103535	Titanium Acetabular Screw 6.5 x 40mm
103536	Titanium Acetabular Screw 6.5 x 45mm
103537	Titanium Acetabular Screw 6.5 x 50mm
103538	Titanium Acetabular Screw 6.5 x 60mm
103539	Titanium Acetabular Screw 6.5 x 70mm

Specialized Instrumentation

"Orthopedic manual surgical instruments are class I exempt per 21 CFR Part 888, section 888.4550.

Part Number	Description
31-103628	M²a™ Trial Liner 41/32mm
31-103629	M²a™ Trial Liner 41/28mm
31-103633	M²a™ Taper Liner Impactor 28mm
31-131004	M²a™ Taper Liner Impactor 32mm
31-103634	M²a™ Liner Extractor 37 Taper
31-103635	M²a™ Liner Extractor 41 Taper
31-434540	M²a™ Ringloc® Cup Inserter
31-434545	Dial-a-Version
31-434543	Dial-a-Version Locking Mechnism





# Sample Labels

REF 10-111148 LOT 123123 2A(TM) / M2A(TM) SYSTEM W/PLUG TI 6AL 4V ALLOY/POROUS COATED WARNING! USE ONLY WITH 37MM O.D. TAPERED LINERS

[LOT] 123123 QTY. 1

BIOMET ORTHOPEDICS, INC. 56 EAST BELL DRIVE
P.O. BOX 587
WARSAW, IN 46581 USA
STERILE R
2004-08
EXPIRY DATE:
2014-08

LOT 123121

REF. 15-105000 M2A(TM) TAPER LINÉR CO-CR-MO ALLOY USE 11-163660/66 M2A(TM) MODULAR HEADS ONLY

**LOT** 123121 QTY. 1 **CE** 0086

BIOMET ORTHOPEDICS, INC. 56 EAST BELL DRIVE P.O. BOX 587 WARSAW, IN 46581 USA STERILE R

2004-08

EXPIRY DA1E: 
2014-08

LOT 123123

11-163660 TYPE I TAPER/METAL ON METAL CO-CR-MO ALLOY



BIOMET ORTHOPEDICS, INC. 56 FAST BEU. DRIVE P.O. BOX 587 WARSAW, IN 46581 USA STERILE R 2004-08

LXPIRY DATE: 
2014-08 ⚠

Biomet Orthopedics, Inc. 56 East Bell Drive P.O. Box 587 Warsaw, Indiana 46581 USA 01-50-0959 Date: 10/04

# Biomet® $M^2a^{TM}$ — Taper Prostheses

#### **Attention Operating Surgeon**

#### DESCRIPTION

The Biomet Metal on Metal Hip Joint Replacement Prosthesis is intended for cemented or uncemented use in primary and revision hip joint replacement procedures. The metal liners are intended for use with specific metal on metal femoral articulating heads. The specialized femoral heads and metal on metal liners are to be used with Biomet primary and revision femoral components.

#### Materials

Femoral Heads CoCrMo Alloy
Acetabular Shells Titanium Alloy
Acetabular Liners CoCrMo Alloy
Porous Coating Titanium Alloy

#### INDICATIONS

- Noninflammatory degenerative joint disease including osteoarthritis and avascular necrosis.
- 2) Rheumatoid arthritis
- 3) Correction of functional deformity.
- Revision procedures where other treatment or devices have been unsuccessful.
- 5) Treatment of non-union, femoral neck fracture, and trochanteric fractures of the proximal femur with head involvement, unmanageable using other techniques.

Patient selection factors to be considered include: 1) need to obtain pain relief and improve function, 2) ability and willingness of the patient to follow instructions, including control of weight and activity level, 3) a good nutritional state of the patient, and 4) the patient must have reached full skeletal maturity.

#### CONTRAINDICATIONS

Absolute contraindications include: infection, sepsis, and ostcomyelitis. Relative contraindications include: 1) uncooperative patient or patient with neurologic disorders who are incapable of following directions, 2) osteoporosis, 3) metabolic disorders which may impair bone formation, 4) osteomalacia, 5) distant foci of infections which may spread to the implant site, 6) rapid joint destruction, marked bone loss or bone resorption apparent on roentgenogram, 7) vascular insufficiency, muscular atrophy, or neuromuscular disease.

#### WARNINGS

Improper selection, placement, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may lead to subsequent reduction in the service life of the prosthetic components. Malalignment of the components or inaccurate implantation can lead to excessive wear and/or failure of the implant or procedure. Inadequate preclosure cleaning (removal of surgical debris) can lead to excessive wear. Improper preoperative or intraoperative implant handling or damage (scratches, dents, etc.) can lead to crevice corrosion, fretting, fatigue fracture and/or excessive wear. Do not modify implants. The surgeon is to be thoroughly familiar with the implants and instruments, prior to performing surgery.

- Use Biomet metal on metal acetabular liners with specified Biomet metal on metal femoral heads.
- Do not use a metallic impactor or punch against the tapered flange of the Titanium shell. Damage to the taper can cause the locking mechanism to malfunction.

- 3) When inserting the final metal liner, do not scratch the outer taper of the metal liner or the inner taper of the metal shell.
- 4) Firmly seat modular head components to prevent dissociation. Thoroughly clean and dry taper prior to attachment of the modular head component to avoid crevice corrosion and improper seating.
- 5) Use only low profile acetabular dome screws with the metal on metal acetabular shells. Acetabular screws are to be fully seated to assure stable fixation and to avoid interference with the metal liner component.
- 6) Prior to seating the liner into the shell component, all surgical debris (tissue fragments, etc.) must be removed from the interior of the shell component, as debris may inhibit the locking mechanism from engaging and securing the liner into the shell component.
- 7) Perforation entirely through the pelvic bone with dome fixation screws or rim screws is to be completely avoided. Caution is to be used when determining and selecting the length of screws to be used, as perforation through the pelvic bone with screws that are too long can cause damage to body structures (blood vessels, etc.) located on the interior side of the pelvis.
- 8) Tight fixation of all non-cemented components at the time of surgery is critical to the success of the procedure. Each component must properly press fit into the host bone which necessitates precise operative technique and the use of specified instruments. Bone stock of adequate quality must be present and appraised at the time of surgery.
- Complete preclosure cleaning and removal of surgical debris at the implant site is critical to minimize wear of the implant articular surfaces.

Biomet joint replacement prostheses provide the surgeon with a means of reducing pain and restoring function for many patients. While these devices are generally successful in attaining these goals they cannot be expected to withstand the activity levels and loads of normal healthy bone and joint tissue.

Accepted practices in postoperative care are important. Failure of the patient to follow postoperative care instructions involving rehabilitation can compromise the success of the procedure. The patient is to be advised of the limitation of the reconstruction and the need for protection of the implants from full load bearing until adequate fixation and healing have occurred. Excessive activity, trauma and weight gain have been implicated with premature failure of the implant by loosening, fracture, and/or wear. Loosening of the implants can result in increased production of wear particles, as well as accelerate damage to bone making successful revision surgery more difficult. The patient is to be made aware and warned of general surgical risks, possible adverse effects as listed, and to follow the instructions of the treating physician including follow-up visits.

#### PRECAUTIONS

Specialized instruments are designed for Biomet joint replacement systems to aid in the accurate implantation of the prosthetic components. The use of instruments or implant components from other systems can result in inaccurate fit, sizing, excessive wear and device failure. Intraoperative fracture or breaking of instruments has been reported. Surgical instruments are subject to wear with normal usage. Instruments, which have experienced extensive use or excessive force, are susceptible to fracture. Surgical instruments should only be used for their intended purpose. Biomet recommends that all instruments be regularly inspected for wear and disfigurement.

Do not reuse implants. While an implant may appear undamaged, previous stress may have created imperfections that would reduce the service life of the implant. Do not treat patients with implants that have been, even momentarily, placed in a different patient.

#### POSSIBLE ADVERSE EFFECTS

Material sensitivity reactions. Implantation of foreign material in tissues
can result in histological reactions involving various sizes of
macrophages and fibroblasts. The clinical significance of this effect is
uncertain, as similar changes may occur as a precursor to or during the
healing process. Particulate wear debris and discoloration from metallic
and polyethylene components of joint implants may be present in
adjacent tissue or fluid. It has been reported that wear debris may

- initiate a cellular response resulting in osteolysis or osteolysis may be a result of loosening of the implant.
- 2) Early or late postoperative, infection, and allergic reaction.
- Intraoperative bone perforation or fracture may occur, particularly in the presence of poor bone stock caused by osteoporosis, bone defects from previous surgery, bone resorption, or while inserting the device.
- 4) Loosening or migration of the implants can occur due to loss of fixation, trauma, malalignment, bone resorption, excessive activity.
- Periarticular calcification or ossification, with or without impediment of joint mobility.
- Inadequate range of motion due to improper selection or positioning of components.
- 7) Undesirable shortening of limb.
- Dislocation and subluxation due to inadequate fixation and improper positioning. Muscle and fibrous tissue laxity can also contribute to these conditions.
- Fatigue fracture of components can occur as a result of loss of fixation, strenuous activity, malalignment, trauma, non-union, or excessive weight.
- 10) Fretting and crevice corrosion can occur at interfaces between components.
- 11) Wear and/or deformation of articulating surfaces.
- 12) Trochanteric avulsion or non-union as a result of excess muscular tension, early weight bearing, or inadequate reattachment.
- 13) Problems of the knee or ankle of the affected limb or contralateral limb aggravated by leg length discrepancy, too much femoral medialization or muscle deficiencies.
- 14) Postoperative bone fracture and pain.
- 15) Elevated metal ion levels have been reported with metal on metal articulating surfaces. Although mechanical testing demonstrates that metal on metal articulating surfaces produce a relatively low amount of particles, the total amount of particulate produced in vivo throughout the service life of the implants remains undetermined. The long-term biological effects of the particulate and metal ions are unknown.

#### STERILITY

Prosthetic components are sterilized by exposure to a minimum dose of 25kGy of gamma radiation. Do not resterilize. Do not use after expiration date.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.


Comments regarding this device can be directed to Attn: Regulatory Dept., Biomet, P.O. Box 587, Warsaw, IN 46581 USA, Fax: 574-372-1683.

Authorized Representative: Biomet U.K., Ltd.


Waterton Industrial Estates, Bridgend, South Wales CF31 3XA, U.K.

CE0086

# Mild-Jilili's metal-on-metal Jarticulation



THREICH TOURS





MIGURE 2


A small starter reamer prepares the acetabulum for subsequent reamers. A hemispherical reamer is utilized to reach the subchondral bone and attain proper sizing.

# Preoperative Planning (Figure 1)


Accurate preoperative planning and acetabular templating are essential for obtaining a successful outcome. Estimate the acetabular size utilizing the RingLoc $^{\circ}$  Mallory-Head $^{\circ}$  or Universal $^{\circ}$  templates along with the appropriate femoral templates in the A/P view. Mark the femoral head center, the neck resection level and the expected femoral stem size on the 14 x 17 A/P radiograph of the femur.

# Incision and Surgical Exposure

The surgical approach is left to the surgeon's discretion. M²a™-Taper instrumentation is compatible with all routine hip exposures (Figure 1).



Progressive hemispherical reaming prepares the acetabulum for the final prosthesis.



A trial metal frame shell gauge is used to determine final shell diameter.


# Acetabular Reaming


Circumferential exposure of the acetabulum should be obtained prior to beginning reaming. It is essential to remove peripheral soft tissue and any remaining acetabular cartilage, to create bleeding bone on the acetabular surface, and to preserve as much of the subchondral bone as possible.

Begin hemispherical reaming of the acetabulum with an acetabular reamer several sizes smaller than the templated diameter of the acetabulum (Figure 2). Progressive reaming is performed until "bleeding" bone, reaming "chatter," or desired template sizing is achieved. In general, the reaming should progress to the floor of the acetabular fossa; under-reaming by Imm may enhance the fit of the final prosthesis, according to the surgeon's judgement of bone stock (Figure 3).

Acetabular trial gauges may be used throughout the reaming process and serve several purposes (Figure 4). The trial gauges determine both reaming accuracy and diameter of the final prosthesis. In addition, the acetabular trials allow judgement of the final position of the implant relative to the peripheral rim. A nerve hook may be inserted into the trial gauge to determine the apposition of the implant to the bone.

During the reaming process, it is important to maintain the structural integrity of the anterior and posterior columns of the acetabular rim. The thickness of the anterior and posterior columns should be monitored by the surgeon to avoid compromise of the bone stock.






# Acetabular Component Insertion

The M²a<sup>™</sup>-Taper system utilizes two different acetabular shells. The Mallory-Head® Radial Shell is a solid porous coated acetabular shell with 8 radial fins that provide excellent rotational stability. The Universal® Shell is a hemispheric implant with two screw holes for additional screw fixation. A 17 degree rim flare provides excellent stability in the rim of the acetabulum.

The acetabular shell component is impacted utilizing the acetabular shell inserter. The shell inserter is threaded directly into the dome of the shell (Figure 5). In order to achieve correct positioning, the Dial-a-Version guide is positioned over the shaft of the acetabular inserter before threading the inserter into the shell and the appropriate anteversion is dialed into the alignment guide. In a lateral position, the inclination rod is perpendicular to the spine to obtain correct horizontal inclination of 45–50 degrees. The anteversion rod is parallel to the spine and should be at 10–15 degrees of anteversion (Figures 6 & 7).

Important: The position of the shell is crucial in the use of the M²a"-Taper Metal-on-Metal Articulation to reduce the risk for impingement.





After the acetabular shell is fully seated, the shell should be checked again to ensure proper orientation. Stability of the bone-implant interface should be checked by applying moderate force to several areas of the rim of the prosthesis. Care must be used to avoid scratching the taper region of the interior of the implant. The acetabular implant should be firmly fixed within the acetabulum, with no gaps between the shell and the acetabulum (Figure 8). If the shell rotates within the acetabulum, a larger shell must be selected and the bone preparation process should be repeated by reaming to the larger size.

Important: If the surgeon desires a change in either the anteversion or inclination, an impactor or punch must not be used against the taper region of the titanium shell. This may damage the taper and cause the connection to malfunction. Avoid contacting the tapered flange with other instrumentation. If the tapered flange is scratched or damaged, the damaged implant should be removed and replaced with a new prosthesis.






#### Acetabular Screw Insertion

Low profile dome screws may be used for adjunctive fixation when utilizing the Universal® 2-Hole Acetabular component and when judged necessary by the surgeon. Screw placement must be chosen carefully to avoid injury to neurovascular structures (Figures 9 & IO). Care should also be exercised when supplement screw fixation is required to avoid damaging or scratching the taper.

Important: Only 6.5mm low profile dome screws may be used with the M<sup>2</sup>a<sup>--</sup>Taper Universal<sup>®</sup> acetabular shells.





### Acetabular Liner Insertion

Following the acetabular shell component insertion and fixation, a trial liner should now be inserted (Figure II). The trial liner is easily inserted into and removed from the metal on metal shell. If preferred, a 3.5mm Hex screwdriver (part number 424496) can be used to tighten down the screw of the trial liner into the apical hole of the acetabular shell. The final liner should be inserted by hand and the surgeon should feel around the face of the shell with his finger to ensure that the liner is aligned properly. Several firm impactions of the acetabular liner are necessary utilizing the liner impactor to ensure stable seating of the device (Figure 12).

Important: When inserting the metal acetabular liner, the interior of the acetabular shell should be carefully cleaned and dried. The taper region of the metal acetabular liner should also be dry before insertion into the acetabular shell. Care should be taken not to scratch the taper surface of the metal liner or the inner taper of the metal shell.



### Modular Head Selection and Impaction

With the acetabular trial or final liner in place, and upon completion of femoral reconstruction, a trial reduction should be performed to confirm restoration of limb length and stability of the hip in all planes.


Based on the chosen final trial head component, the corresponding modular head must now be selected. Only modular heads labeled for the  $M^2a^{TM}$ -Taper articulation should be utilized. The precision tolerances of these femoral heads are necessary for optimum wear (Figure 13).

# Important: Be sure taper surfaces are clean and dry before seating the modular head on the stem taper.

Impact the modular head onto the stem with several brisk mallet strikes using a plastic head impactor only. Metal impactors or any other metallic objects may scratch the modular head bearing surface or flatten out the radius in this region and, therefore, should not be used. If the modular head becomes scratched, it must be replaced.



Above photo shows the retraction of the snap cylinder used to disassociate liner from shell.



# Separation of the Taper Junction

If the Metal-on-Metal Liner needs to be removed from the Metal-on-Metal Taper Acetabular Shell, the Metal-on-Metal Taper Liner Extractor is to be used. Turn the handle of the liner extractor to the left to expand the prongs and align with two opposing grooves on the acetabular shell (Figure 14). To ensure proper extension of the prongs, firmly tap the butt of the handle. Once the liner extractor is in contact with the shell, turn the handle to the right tightening the prongs around the liner. Once the outer hub is tightened against the acetabular shell, pull the snap cylinder on the liner extractor and release (Figure 14)\* Remove the liner from the shell and turn the handle to the left to release the liner from the extractor (Figure 15).

\*The vibration created from the impact of the snap cylinder will loosen the liner from the shell.

# **Predicate Device Comparison Table**

	New device	Metal on Metal Acetabular System	M²a™ 32mm Taper System	Mallory/Head PF
				Acetabular Component
Manufacturer	Biomet	Biomet	Biomet	Biomet
510(k) number	New	K993438	K003363	K861114
Indications for use	1) Non-inflammatory	1) Non-inflammatory	1) Non-inflammatory	None specified
	degenerative joint disease	degenerative joint disease	degenerative joint disease	
	including osteoarthritis and	including avascular necrosis,	including osteoarthritis and	
	avascular necrosis	diastrropic varient, fracture of	avascular necrosis	
	2) Rheumatoid arthritis	the pelvis, fused hip, leg	<ol><li>Rheumatoid arthritis</li></ol>	
	3) Correction of functional	perthes, osteoarthritis, slipped	<ol><li>Correction of functional</li></ol>	
	deformity	capital epiphysis, subcapital	deformity	
	4) Revision procedures where	fractures, and traumatic	4) Revision procedures where	
	other treatment or devices	arthritis	other treatment or devices	
	have been unsuccessful,	2) Rheumatoid arthritis	have failed	
	5) Treatment of non-union,	3) Correction of functional	<ol><li>Treatment of non-union,</li></ol>	
	femoral neck fracture,	deformity	femoral neck fracture,	
	trochanteric fractures of the	4) Treatment of non-union,	trochanteric fractures of the	
	proximal femur with head	femoral neck fracture,	proximal femur with head	
	involvement, unmanageable	trochanteric fractures of the	involvement, unmanageable	
	using other techniques	proximal femur with head	using other techniques	
		involvement, unmanageable		
		using other techniques		
		5) Revision of previously failed		
		total hip arthroplasty		
Intended use	Cemented and Uncemented	Uncemented	Uncemented	Cemented

	New device	Metal on Metal Acetabular System	setabular System	M2a 32mm Taper System	aper System	Mallory/Head PF
						Acetabular Component
510(k) number	New	K993438	1438	K003363	363	K861114
	(1) (1) (1) (1) (1) (1) (1) (1) (1) (1)	Acetabular Liner	ar Liner			
Liner Material	Co-Cr-Mo	Co-Cr-Mo	r-Mo	Co-Cr-Mo	r-Mo	Not Applicable
Liner Internal Diameter	28mm & 32mm	28mm	mn	32mm	nm	
Liner Size Designation	37mm & 41mm	37mm	uu	41mm	nm	
Liner Locking Mechanism	Taper	Taper	ber .	Taper	oer .	
	(1) (1) (1) (1) (1) (1) (1) (1) (1) (1)	Acetabular Shell	lar Shell	A Comment of the Comm		
Shell Style	Mailory/Head	Mallory/Head	Universal	Mallory/Head	Universal	Mallory/Head
Shell Material	Ti-6Al-4V	Ti-6Al-4V	Ti-6Al-4V	Ti-6Al-4V	Ti-6Al-4V	Ti-6Al-4V
Liner Locking Mechanism	Taper	Taper	Taper	Taper	Taper	Hex-Loc
Shell Outer Diameters	48mm - 70mm by 2mm	52mm - 70mm	48mm - 70mm	52mm - 70mm	52mm - 70mm	46mm-70mm by 4mm
	increments	by 2mm	by 2mm	by 2mm	by 2mm	increments
		increments	increments	increments	increments	
Shell Features	Radial Fins	Radial Fins	2 Holes	Radial Fins	2 Holes	Fins
	2 Holes	Solid	Plasma Spray	Solid	Plasma Spray	Plasma Spray
	Plasma Spray	Plasma Spray	Straight Cutouts	Plasma Spray	Straight Cutouts	
	Angled Cutouts	Straight Cutouts		Straight Cutouts		



Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

#### MAY 1 8 2000

Ms. Michelle L. McKinley Regulatory Specialist Biomet, Inc. Airport Industrial Park P.O. Box 587 Warsaw, Indiana 46581-0587

Re: K993438/S1

Trade Name: Metal on Metal Acetabular Component

Regulatory Class: III Product Code: KWA Dated: February 18, 2000 Received: February 22, 2000

Dear Ms. McKinley:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general control provisions of the Act. The general control provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

#### Page 2 - Ms. Michelle L. McKinley

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its Internet address "http://www.fda.gov/cdrh/dsmamain.html".

Sincerely yours,

Celia M. Whitten, Ph.D., M.D.

Through Van

Director

Division of General, Restorative and

Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

	Page of	
510(k) Number if Known: K993438  Device Name: Metal on Metal Acetabular System		
The Metal on Metal Acetabular System is indicated for used replacement due to the following:	in patients requiring total hip	
<ul> <li>a.) Non-inflammatory degenerative joint disease including avascular necrosis, diastrophic variant, fracture of the pelvis, fused hip, leg perthes, osteoarthritis, slipped capital epiphysis, subcapital fractures, and traumatic arthritis</li> <li>b.) Rheumatoid arthritis</li> <li>c.) Correction of functional deformity</li> <li>d.) Treatment of non-union, femoral neck fracture, and trochanteric fractures of the proximal femur with head involvement, unmanageable using other techniques.</li> <li>e.) Revision of previously failed total hip arthroplasty.</li> </ul>		
,		
(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)		
Concurrence of CDRH, Office of Device Evaluation (ODE)		
1:00011ption 000	Over the Counter Use(Optional Format 1-2-96)	
(Division Sign-Off) Division of General Restorative Devices 510(k) Number K 5 5 3 4 3 8		



DEC 1 5 2000

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Ms. Michelle L. Mckinley Regulatory Specialist Biomet, Inc. P.O.Box 587 Warsaw, Indiana 46582

Re: K003363

Trade Name: M2A 32 MM Taper System

Regulatory Class: III Product Code: KWA Dated: December 6, 2000 Received: December 7, 2000

Dear Ms. Mckinley:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general control provisions of the Act. The general control provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

### Page 2 - Ms. Michelle L. Mckinley

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its Internet address "http://www.fda.gov/cdrh/dsmamain.html".

Sincerely yours,

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative and

Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

**Device Name:** M2a<sup>™</sup> 32mm Taper System

### Indications for Use:

The M2a<sup>TM</sup> 32mm Taper System is indicated for use in patients requiring total hip replacement due to the following:

1) Non-inflammatory degenerative joint disease including osteoarthritis and avascular necrosis; 2) rheumatoid arthritis; 3) correction of functional deformity; 4) revision procedures where other treatment or devices have failed; 5) treatment of non-union, femoral neck fracture, and trochanteric fractures of the proximal femur with head involvement, unmanageable using other techniques.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE OF ANOTHER PAGE IS NEEDED)

(Division Sign-Off)

Division of General Restorative Devices K 003363

510(k) Number \_

000007

Ms. Judith Dermody Clinical Research Assistant Biomet Inc. Corporate Headquarters Airport Industrial Park P.O. Box 587 Warsaw, Indiana 46580

APR 17 1986

Food and Drug Administration 8757 Georgia Avenue Silver Spring MD 20910

Re: K861114

Mallory/Head - PF Acetabular

Component

Dated: March 17, 1986 Received: March 25, 1986

Dear Ms. Dermody:

We have reviewed your Section 510(k) notification of intent to market the above device and we have determined the device to be substantially equivalent to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments. You may, therefore, market your device subject to the general controls provisions of the Federal Food, Drug, and Cosmetic Act (Act) until such time as your device has been classified under Section 513. At that time, if your device is classified into either class II (Performance Standards) or class III (Premarket Approval), it would be subject to additional controls. Please note: This action does not affect any obligation you might have under the Radiation Control for Health and Safety Act of 1968, or other Federal Laws or regulations.

General controls presently include regulations on annual registration, listing of devices, good manufacturing practice, labeling, and the misbranding and adulteration provisions of the Act. In the future, the scope of general controls may be broadened to include additional regulations.

All regulations and information on meetings of the device advisory committees, their recommendations, and the final decisions of the Food and Drug Administration (FDA) will be published in the Federal Register. We suggest you subscribe to this publication so you can convey your views to FDA if you desire and be notified of any additional requirements imposed on your device. Subscriptions may be obtained from the Superintendent of Documents, U.S. Government Printing Office, Washington, D.C. 20402. Such information also may be reviewed in the Dockets Management Branch (HFA-305), Food and Drug Administration, Room 4-62, 5600 Fishers Lane, Rockville, Maryland 20857.

This letter does not in any way denote official FDA approval of your device or its labeling. Any representation that creates an impression of official approval of this device because of compliance with the premarket notification regulations is misleading and constitutes misbranding. If you desire advice on the labeling for your device or other information on your responsibilities under the Act, please contact the Office of Compliance, Division of Compliance Operations (HFZ-320), 8757 Georgia Avenue, Silver Spring, Maryland 20910.

Sincerely yours,

Cal A. Foron

Carl A. Larson, Ph.D.
Director
Division of Surgical
and Rehabilitation Devices
Office of Device Evaluation
Center for Devices
and Radiological Health

### Engineering Justification for $M^2 a^{TM}/C^2 a^{TM}$ Shells


The  $M^2a^{TM}/C^2a^{TM}$  Shells are substantially equivalent to our currently cleared  $M^2a^{TM}$  Shells (K003363). Both series of shells have the same dimensions for the locking cone interface. The size range is also equivalent in both series of shells. Both series of shells also accept the same Co-Cr-Mo metal bearing surface. The exterior geometry of the shells is identical.


The dimensions for the locking mechanism on the M²aTM/C²aTM shells are both an 18°55'00" taper. The taper engagement length is also identical to the M²aTM Shells. The only difference in the locking mechanism is the surface finish of the female tapers. The M²aTM system has a smoother female taper surface finish, whereas the M²aTM/C²aTM shells have a roughened female taper surface because the system is designed for future use with Biolox® Forte Alumina Inserts from CeramTec AG as well as the Co-Cr-Mo Inserts. The surface finish needed to be roughened in order to distribute the stress across the insert. This rougher surface acts like small peaks that can be conformed to the surface of the mating component actually increasing the locking strength of the interface as well as distributing the stress over more surface area. Ceramic Inserts were tested in a push-out (MT3045 and MT-3260) scenario with the smoother surface as well as the roughened surface. The smoother surface yielded an average push-out force of 180 lbs. The roughened surface yielded an average push-out force of 180 lbs. The roughened surface yielded an average push-out force of 180 lbs. Both of these values far exceed the requirements set by CeramTec AG for ceramic liner push-out which is only 45lbs.


Sizing for the  $M^2a^{TM}/C^2a^{TM}$  shell system is 48mm to 70mm in increments of 2mm this is identical to the cleared  $M^2a^{TM}$  system. The exterior geometry of the system is based off of our Mallory/Head Radial® design with 8 evenly spaced fins. This geometry is also used in the  $M^2a^{TM}$  system. The metal insert the shells accept is also identical to the  $M^2a^{TM}$  System. The 48mm and 50mm shells accept a 28/37 tapered insert and the 54mm-70mm shells accept a 32/41 tapered insert.


Based on the previous comparison between the two shell systems I believe the  $M^2 a^{TM}/C^2 a^{TM}$  shell system is substantially equivalent to the  $M^2 a^{TM}$  system and no further testing is required.


W. Vasor Stone

















Conclusions	
(b)(4)	
References	
(b)(4)	


# Appendix 1

Raw Data













## Appendix 2


Biomet Laboratory Request





















## Appendix 3


Drawings

























### Results









## References





# Appendix 1


Raw Data







































# Appendix 2

Biome Laboratory Request












# Appendix 3


Drawings


118
















### 510(k) Summary

**Applicant/Sponsor:**Biomet Manufacturing Corp.

Contact Person: Patricia Sandborn Beres

Senior Regulatory Specialist

**Proprietary Name**: M²a™/C²a™ Acetabular System

Common Name: Metallic Acetabular System

### **Classification Name:**

1. Hip joint metal/metal semi-constrained, with an uncemented acetabular component prosthesis (21 CFR 888.3330)

2. Hip joint metal/metal semi-constrained, with a cemented acetabular component prosthesis (21 CFR 888.3320)

Legally Marketed Devices To Which Substantial Equivalence Is Claimed: Biomet devices:

K993438 - Metal on Metal Acetabular System

K003363 - M²a™ 32mm Taper System

K861114 - Mallory/Head PF Acetabular Component

**Device Description:** The M²a™/C²a™ Acetabular System consists of a titanium outer acetabular shell with a cobalt alloy metallic liner for metal on metal articulation.

The acetabular shells are hemispherical in shape to closely match the natural acetabulum. Two screw holes in the dome allow for additional fixation by the use of screws. The outer surface of the shells are covered with Biomet's plasma sprayed coating.

The metallic cobalt alloy bearing liner fits into the outer shell by means of a taper similar to the taper used for the attachment of a modular head to a femoral stem. The metallic liners articulate with cobalt alloy modular heads.

MAILING ADDRESS P.O. Box 587 Warsaw, IN 46581-0587 SHIPPING ADDRESS 56 E. Bell Drive Warsaw, IN 46582 M²a™/C²a™ Acetabular System 510(k) Summary Page 2

**Intended Use:** The M²a™/C²a™ Acetabular System is intended for cemented or noncemented use in cases of:

- 1) Non-inflammatory degenerative joint disease including osteoarthritis and avascular necrosis
- 2) Rheumatoid arthritis
- 3) Correction of functional deformity
- 4) Revision procedures where other treatment or devices have been unsuccessful
- 5) Treatment of non-union, femoral neck fracture, trochanteric fractures of the proximal femur with head involvement, unmanageable using other techniques

**Summary of Technologies:** The technological characteristics of the new device are similar of identical to the predicates.

Non-Clinical Testing: None provided

Clinical Testing: None provided.

All trademarks are property of Biomet, Inc.

Pages 127 through 243 are not included as part of this submission

M<sup>2</sup>a<sup>™</sup>/C<sup>2</sup>a<sup>™</sup>Acetabular System 510(k) Premarket Notification Biomet Manufacturing Corp.
October 13, 2004

Appendix 1 & 2

TBA/CONFORM D 1.30

### **Class III Certification and Summary**

[As required by 21 CFR 807.94]

I certify in my capacity as Senior Regulatory Specialist of Biomet Manufacturing Corp, that I have conducted a reasonable search of all information known or otherwise available about the types and causes of safety and effectiveness problems that have been reported for the metal on metal acetabular components. I further certify that I am aware of the types of problems to which the metal on metal acetabular components are susceptible and that, to the best of my knowledge, the following summary of the types and causes of safety and effectiveness problems about the metal on metal acetabular components is complete and accurate.

Patricia 5 Beur Signature

Patricia Sandborn Beres\_\_\_\_

Typed Name

10/6/04 Date

### **Class III Device Summary**

### Metal on Metal Articulating Surfaces

I. Medical Device Reports/Vigilance Reports: Summary of MAUDE Medical Device Reports

A reasonable effort was made to find all adverse reports made for these devices under the Medical Device Reporting (MDR) regulations and under the vigilance-reporting requirement for medical devices under Article 10 of the European Medical Device Directive (MDD). A search of publicly available information yielded thirty-three (33) reports filed with the FDA for the metal/metal semi-constrained total hip prostheses. A summary of the reports is located on the next page, followed by MAUDE database reports of each incident listed.

# SUMMARY OF ADVERSE EVENTS TABLE

Adverse	Disassociation of anti-rotational pin from insert.	Deep infection resulting in explant of device.	Revision due to impingement between femoral stem and acetabular insert.	Pin insert came out after 1.5 years.	Multiple dislocations, patient revised 3 times.	Disassociation of Metasul insert from APRII shell. Patient also experienced two heavy falls previously.	Doctor impacted the insert into the shell but it did not seat. When he pulled the insert out to inspect, the locking pin was missing and located in one of the slots. Resulted in a 30-minute delay in surgery.	Patient complaining of hearing a "pop" and feeling a "pop" and thinks they have dislocated their hip. Physician notified.	Revision of femoral head at 7 weeks.	Metasul insert dislocated after 4 years implantation, significant metalosis. Severe pain.	APR Metasul insert disengaged nine months after surgery.	Disassociation of liner from cup.	Anti-rotational pin dislocated from the insert.	Surgery time extended due to failure to properly assemble acetabular liner into shell.	Doctor reported patient with pain, disassociation of liner from shell.	Disassociation of APR Metasul acetabular insert, size 59mm. Patient fell 2001.	APR Metasul acetabular liner disassociation, patient revised in 2002.
Date FDA Received Notification	5/24/2000	7/5/2000	8/17/2000	12/1/2000	1/5/2001	2/22/2001	6/8/2001	6/13/2001	7/6/2001	7/18/2001	9/18/2001	12/18/2001	1/7/2002	1/9/2002	2/6/2002	2/27/2002	3/19/2002
Date of Report	4/24/2000	6/16/2000	7/21/2000	11/6/2000	12/8/2000	1/23/2001	5/10/2001	5/14/2001	6/8/2001	6/18/2001	8/23/2001	11/26/2001	12/12/2001	12/17/2001	1/28/2002	1/28/2002	2/27/2002
Device 510K Number	K993569	K674728	K993569	K993569	K993569	K993569	K993569	K974728	K974728	K993569	K993269	K993569	K993569	K003250	K993569	K993569	K993569
Device Manufacturer	Sulzer Orthopedic, Inc	Sulzer Orthopedic, Inc	Sulzer Orthopedic, Inc	Sulzer Orthopedic, Inc	Sulzer Orthopedic, Inc	Sulzer Orthopedic, Inc	Sulzer Orthopedic, Inc	Sulzer Orthopedic, Inc	Sulzer Orthopedic, Inc	Sulzer Orthopedic, Inc	Sulzer Orthopedic, Inc	Sulzer Orthopedic, Inc	Sulzer Orthopedic, Inc	Encore Orthopedics, Inc	Sulzer Orthopedic, Inc	Sulzer Orthopedic, Inc	Sulzer Orthopedic, Inc
Event:	٦	2	3	4	5	9	2	&	6	10	7	12	13	14	15	16	17

. 'Adverse Event	Disassociation of APR Metasul insert.	Reported pain and breakage of implant.	Disassociation of insert 2 years 9 months after initial surgery.	Patient was revised due to dislocation of the hip.	Eleven total joint and/or revision procedures performed with five of these cases resulting in postoperative infection. Various implant components implanted. Concluded not to be related to implants because similar reports of infections have not been reported from other user facilities for any of the listed manufacturing lots.	Second revision required because of instability of the hip.  First revision was on a recall shell.	Following total hip arthroplasty performed in 2002, patient continued to experience hip pain. Patient underwent additional surgery 1 month later, where osteophytes were removed and modular head component was exchanged.	The package was split – inner plastic portion.	The metal insert was locked into the shell off-axis. Insert could not be removed from acetabular shell.	Patient felt pain due to dislocation. Patient was revised.	Dislocation of APR Metasul insert.	Patient had a dislocation and was revised.	Delay in surgery.	Patient was revised.	Patient was revised in 2003.	Metal insert did not seat properly. Surgery was extended by 30 minutes due to the product problem.	Ring broken, worn liner. Revised liner only to non- constrained.
Date FDA E Received	4/12/2002	4/12/2002	5/1/2002	6/26/2002	11/20/2002	12/17/2002	1/7/2003	1/10/2003	3/27/2003	6/25/2003	7/16/2003	7/17/2003	10/14/2003	12/16/2003	12/23/2003	1/16/2004	1/09/2004
Date of Report	3/15/2002	3/15/2002	4/15/2002	6/4/2002	11/8/2002	12/10/2002	12/9/2002	1/10/2003	3/27/2003	6/10/2003	6/19/2003	6/25/2003	10/14/2003	12/16/2003	12/17/2003	12/19/2003	12/10/2003
** Device **  *** 510K**  Number	3	K993569	K993569	K974728	K011110	K974728	K011110	K003523	K003523	K993569	K993569	K993569	K974728	K993569	K993569	K003523	K924492
P. Device Manufacturer	Sulzer Orthopedic, Inc	Sulzer Orthopedic, Inc	Sulzer Orthopedic, Inc	Sulzer Orthopedic, Inc	Biomet Inc.	Centerpulse Orthopedics, Inc	Biomet, Inc	DePuy International, LTD	DePuy International, LTD	Centerpulse Orthopedics, Inc	Centerpulse Orthopedics, Inc	Centerpulse Orthopedics, Inc	Centerpulse Orthopedics, Inc	Centerpulse Orthopedics, Inc	Centerpulse Orthopedics, Inc	DePuy Orthopaedics, Inc	DePuy Orthopedics
Event Sequence	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34

Adverse Event	Attorney reports this patient was revised due to premature polyethylene failure.	Revision	Surgeon claims polyethylene wear started to appear on x-rays 3 years after THA	Patient was revised in 2004
Date FDA Received Notification	2/24/2004	2/28/2003	5/25/2004	6/24/2004
Date of Report	2/10/2004	2/13/2003	5/19/2004	6/24/2004
Device 510K	K924492	K003758	K924492	K993569
Device Manufacturer	DePuy-Raynham	Centerpulse Orthopedics Ltd.	DePuy Orthopedics, Inc.	Zimmer Austin, Inc.
Events. Sequence	35	36	37	38

## U.S. Food and Drug Administration CENTER FOR DEVICES AND RADIOLOGICAL HEALTH



FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 1

New Search

Back To Search Results

## Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X55 METASUL APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X55 METASUL APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-055

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

SULZER ORTHOPEDICS INC.

Manufacturer (Section F) 9900 Spectrum

Austin TX 78717

SULZER ORTHOPEDICS INC.

Manufacturer (Section D) 9900 Spectrum

Austin TX 78717

Randy Jasek, Supervisor

9900 Spectrum Dr

Manufacturer Contact Austin , TX 78717

(512) 432 -9611

Device Event Key 269992

MDR Report Key 278978

**Event Key** 261599

Report Number 2935620-2000-00012

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Event Type Injury

Type of Report Initial, Followup

Report Date 04/24/2000

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 05/24/2000

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-055

Device LOT Number 1251199

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Date Returned to Manufacturer 04/24/2000

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 04/24/2000

Was Device Evaluated By Manufacturer? Yes

Date Device Manufactured 08/01/1996

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization

Adverse Event or Product Problem Description

Report Date: 04/24/2000 MDR Text Key: 953779 Patient Sequence Number: 1
Allegedly the anti-rotation pin became dislodged from the polyethylene acetabular insert.

Additional Manufacturer Narrative

Report Date: 04/24/2000 MDR Text Key: 1012481

H6 method reviewed mfg/inspection records with no form, fit or function discrepancies noted. H6 conclusions the cause of dissociation of the anti-rotation pin from the insert cannot be determined due to insufficient info provided.

#### Patient TREATMENT DATA

Date Received: 08/08/2000 Patient Sequence Number: 1

# Treatment Treatment Date

1 3442 28MM +4MM ALLOPRO BATORY FEM HD

01/01/1998

2 (LOT# B069262)(10/98).

01/01/1998

Database contains data received through March 30, 2004

Center for Devices and Radiological Health / CDRH

## U.S. Food and Drug Administration GENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 2

New Search

Back To Search Results

## Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name SZ 28X57MM STD METASUL INS I-O

Type of Device HIP PROSTHESIS

Baseline Brand Name SZ 28X57MM STD METASUL INS I-O

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4372-28-057

Baseline Device Family NA

Baseline Device 510(K) Number K974728

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 08/03/1999

SULZER ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr. Austin TX 78717

Randy Jasek, Supervisor

Manufacturer Contact

9900 Spectrum Dr Austin , TX 78717

(512) 432 -9611

Device Event Key 275284

MDR Report Key 284467

**Event Key** 266856

Report Number 2935620-2000-00022

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Event Type Injury

Type of Report Initial, Followup

Report Date 06/16/2000

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 07/05/2000

is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4372-28-057

Device LOT Number 1330653

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Date Returned to Manufacturer 06/19/2000

is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 06/16/2000

Was Device Evaluated By Manufacturer? Yes

Date Device Manufactured 10/01/1998

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Required Intervention

### Adverse Event or Product Problem Description

Report Date: 06/16/2000 MDR Text Key: 974574 Patient Sequence Number: 1

Pt weight bearing for only one week only then developed debilitating pain and returned to bed until deep infection was determined and explanted in 2000.

#### Additional Manufacturer Narrative

Report Date: 06/16/2000 MDR Text Key: 1026547

H6 method reviewed (other) reviewed mfr/inspection records, with no discrepancies noted. H6 conclusions (other) the revision surgery was due to deep septic infection, and was unrelated to sulzer orthopedics inc. Implants. H10 documentation provided from agent states "revision surgery was due to deep septic infection and was unrelated to zulzer orthopedics inc. Implants. ".

#### Patient TREATMENT DATA

Date Received: 08/24/2000 Patient Sequence Number: 1

#	Treatment	Treatment Date
1 7356-01-102 NATURAL HIP HA COLLARLESS SZ 2 LT		11/01/1999
2 (LOT # 1370384)(11/1999)		11/01/1999
3 7340-28-400 METASUL COCR HD 12/14 +4 NECK 28MM		11/01/1999
4 (LOT # 1340044)(11/1999).		11/01/1999

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

## U.S. Food and Drug Administration GENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 3

New Search

Back To Search Results

### Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X53 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X53 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-053

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

SULZER ORTHOPEDICS INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Randy Jasek, Supervisor

9900 Spectrum Dr Manufacturer Contact

Austin, TX 78717

(512) 432 -9611

Device Event Key 281298

MDR Report Key 290650

Event Key 272806

Report Number 2935620-2000-00030

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Event Type Injury

Type of Report Initial, Followup

Report Date 07/21/2000

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 08/17/2000

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-053

Device LOT Number 1187760

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Date Returned to Manufacturer 08/09/2000

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 07/21/2000

Was Device Evaluated By Manufacturer? No

Date Device Manufactured 09/01/1995

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization Required Intervention

#### Adverse Event or Product Problem Description

Report Date: 07/21/2000 MDR Text Key: 997321 Patient Sequence Number: 1

It was reported: revision hip surgery was performed due to impingement between the femoral stem and the acetabular insert.

#### Additional Manufacturer Narrative

Report Date: 07/21/2000 MDR Text Key: 1064198

Added section h6. Corrected section h3. H6: method: review of mfr/inspection records with no discrepancies noted. H6: results: indications are that early impingement was occurring, possibly with soft tissue or improperly positioned acetabular or femoral components. H6: results: possible surgical technique. H6: conclusions: the cause of occurrence could not be determined since x-ray and the femoral stem were not provided. Speculation is that the femoral stem impinged on the acetabular insert.

#### Patient TREATMENT DATA

Date Received: 10/05/2000 Patient Sequence Number: 1

#		Treatment	Treatment Date
1	1. 4310-02-053 APR II 12 SLOT SHELL 53MM (LOT#		01/01/2000
2	1195949)(11/2000).		01/01/2000
3	2. 7340-28-004 METASUL COCR HD 12/14-4 NECK 28MM		01/01/2000
4	(LOT#1181955)(11/2000)		01/01/2000
	256		

5 3. 7354-01-203 NATURAL-HIP POR COLL STM LT SZ 3

6 (LOT# 1213577)(11/2000).

01/01/2000 01/01/2000

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer EDA Home Page | Search EDA Site | EDA A-Z Index | Contact FDA | HHS Home Page

## U.S. Food and Drug Administration CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 4

New Search

Back To Search Results

### Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X49 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X49 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-049

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

SULZER ORTHOPEDICS, LTD.

Grabenstrasse 25 Manufacturer (Section F)

Baar

SWITZERLAND CH-6341

SULZER ORTHOPEDICS, LTD.

Grabenstrasse 25 Manufacturer (Section D)

Baar

SWITZERLAND CH-6341

Randy Jasek, Supervisor

9900 Spectrum Dr Manufacturer Contact

Austin , TX 78717 (512) 432 -9611

Device Event Key 296783

MDR Report Key 306697

Event Key 288203

Report Number 2935620-2000-00062

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Event Type Injury

Type of Report Initial, Followup

Report Date 11/06/2000

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 12/01/2000

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-049

Device LOT Number 1303668

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Date Returned to Manufacturer 11/13/2000

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 11/06/2000

Was Device Evaluated By Manufacturer? Yes

Date Device Manufactured 09/01/1997

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Required Intervention

Adverse Event or Product Problem Description

Report Date: 11/06/2000 MDR Text Key: 1055495 Patient Sequence Number: 1

It was reported: the pin in the metasul insert came out after 1. 5 years.

Additional Manufacturer Narrative

Report Date: 11/06/2000 MDR Text Key: 1156562

H6: method: reviewed mfr/inspection records, with no discrepancies noted. H6: conclusions: insufficient info was provided to determine the actual cause of this complaint.

Patient TREATMENT DATA

Date Received: 04/25/2001 Patient Sequence Number: 1

Treatment

**Treatment Date** 

1 3461 AP BALL HEAD 28/14 (LOT# B285607)(2000).

01/01/2000

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

## U.S. Food and Drug Administration CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 5

New Search

Back To Search Results

## Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X49 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X49 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-049

Baseline Device Family NA

Baseline Device 510(K) Number K993569

is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

SULZER ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Randy Jasek, Supervisor

9900 Spectrum Dr Manufacturer Contact

Austin, TX 78717

(512) 432 -9611

Device Event Key 300804

MDR Report Key 310957

**Event Key** 292272

Report Number 2935620-2000-00075

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Event Type Injury

Type of Report Initial

Report Date 12/08/2000

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 01/05/2001

is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-049

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Date Returned to Manufacturer 12/08/2000

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 12/08/2000

Was Device Evaluated By Manufacturer? No

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Required Intervention

**Adverse Event or Product Problem Description** 

Report Date: 12/08/2000 MDR Text Key: 1071225 Patient Sequence Number: 1

It was reported: pt underwent total hip arthroplasty (tha) in 1998. Subsequently the pt was revised 3 times due to dislocations. Pt underwent the last tha in 2000 where the insert and ball head were replaced.

Patient TREATMENT DATA

#

Date Received: 01/05/2001 Patient Sequence Number: 1

Treatment

**Treatment Date** 

1,3462 METASUL HEAD 29L/+4 12/14 (LOT#UNK) (11/2000),

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

6

# U.S. Food and Drug Administration & Health Human

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards

New Search

Back To Search Results

## Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X49 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X49 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-049

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

SULZER ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Randy Jasek, Supervisor

Manufacturer Contact 9900 Spectrum Dr Austin, TX 78717

(512) 432 -9611

Device Event Key 306464

MDR Report Key 316925

**Event Key** 297939

Report Number 2935620-2001-00003

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

169

Event Type Injury

Type of Report Initial

Report Date 01/23/2001

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 02/22/2001

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-049

Device LOT Number 1230114

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Date Returned to Manufacturer 01/25/2001

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 01/23/2001

Was Device Evaluated By Manufacturer? No

Date Device Manufactured 03/01/1996

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization Required Intervention

### Adverse Event or Product Problem Description

Report Date: 01/23/2001 MDR Text Key: 1092275 Patient Sequence Number: 1

It was reported: 11/13/2000: sudden "clunk" in hip-could not walk, x-ray in emergency dept: disassociation of metasul insert from the apr ii shell. Pt experienced two heavy falls, one onto their back in march 2000 and another fall forward in july 2000.

#### Patient TREATMENT DATA

Date Received: 02/22/2001 Patient Sequence Number: 1

Treatment Treatment Date

1 19.28.06 METASUL HEAD 28MM (11/2000)

11/01/2000

2 (LOT# B154629).

11/01/2000

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer EDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

## U.S. Food and Drug Administration

### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 7

New Search

Back To Search Results

## Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X49 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X49 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-049

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

SULZER ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Randy Jasek, Qa Supervisor

9900 Spectrum Drive

Manufacturer Contact Austin, TX 78717

(512) 432 -9611

Device Event Key 325690

MDR Report Key 336384

**Event Key** 316615

Report Number 2935620-2001-01058

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Event Type Injury

Type of Report Initial

Report Date 05/10/2001

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 06/08/2001

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-049

Device LOT Number 1427124

Was Device Available For Evaluation? No

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 05/10/2001

Was Device Evaluated By Manufacturer? Device Not Returned To Manufacturer

Date Device Manufactured 07/01/2000

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization

Adverse Event or Product Problem Description

Report Date: 05/10/2001 MDR Text Key: 1162660 Patient Sequence Number: 1

It was reported: the dr impacted the insert into the shell, but it did not seat. When he pulled the insert out to inspect, staff noticed the locking pin was missing and now in one of the slots. Resulted in a 30 minute delay in surgery. The patient was not at risk.

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search EDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

## U.S. Food and Drug Administration CENTER FOR DEVICES AND RADIOLOGICAL HEALTH



FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 8

New Search

Back To Search Results

## Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name SZ 28X55MM STD METASUL INS I-O

Type of Device HIP PROSTHESIS

Baseline Brand Name SZ 28 X 55MM STD METASUL INS I-O

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4372-28-055

Baseline Device Family NA

Baseline Device 510(K) Number K974728

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 08/03/1999

SULZER ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Randy Jasek, Supervisor

9900 Spectrum Drive

Manufacturer Contact Austin , TX 78717

(512) 432 -9611

Device Event Key 326311

MDR Report Key 337007

Event Key 317178

Report Number 2935620-2001-01060

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

4/26/2004

Source Type Health Professional

Event Type Injury

Type of Report Initial

Report Date 04/27/2001,05/14/2001

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 06/13/2001

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4372-28-055

Device LOT Number 1417750

Was Device Available For Evaluation? No

Is The Reporter A Health Professional? Yes

Was the Report Sent to FDA? No

Distributor Facility Aware Date 04/25/2001

Device Age 1 mo

Event Location Hospital

Date Report TO Manufacturer 04/27/2001

Date Manufacturer Received 05/14/2001

Was Device Evaluated By Manufacturer? Device Not Returned To Manufacturer

Date Device Manufactured 07/01/2000

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization Required Intervention

Adverse Event or Product Problem Description

Report Date: 04/27/2001 MDR Text Key: 1165142 Patient Sequence Number: 1

It was reported: pt has a left total hip replacement. After hanging up the phone pt complaining of hearing a "pop" and feeling a "pop". Pt stated "i think i dislocated my hip again. " physician was notified.

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

## FDA U.S. Food and Drug Administration & Department of Health and Human Services

### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards

9

New Search

Back To Search Results

## Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name COCR HD 12/14 +4MM NECK 28MM METASUL

Type of Device HIP PROSTHESIS

Baseline Brand Name COCR HD 12/14 +4MM NECK 28MM METASUL

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 7340-28-400

Baseline Device Family NA

Baseline Device 510(K) Number K974728

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 08/03/1999

SULZER ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Randy Jasek, Supervisor

Manufacturer Contact 9900 Spectrum Drive

Austin , TX 78717

(512) 432 -9611

Device Event Key 329938

MDR Report Key 340619

Event Key 320719

Report Number 2935620-2001-01186

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

269

Event Type Injury

Type of Report Initial

Report Date 06/08/2001

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 07/06/2001

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 7340-28-400

Device LOT Number 1327368-F

Was Device Available For Evaluation? No

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 06/08/2001

Was Device Evaluated By Manufacturer? Device Not Returned To Manufacturer

Date Device Manufactured 10/01/1998

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization

Adverse Event or Product Problem Description

Report Date: 06/08/2001 MDR Text Key: 1178778 Patient Sequence Number: 1

It was reported: revision of femoral head 7 weeks later.

Database contains data received through March 30, 2004

CDRH, Home Page | CDRH, A.-Z. Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z. Index | Contact FDA | HHS, Home Page

## U.S. Food and Drug Administration

#### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 10

New Search

Back To Search Results

### Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X49 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X49 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-049

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

SULZER ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Randy Jasek, Supervisor

9900 Spectrum Drive Manufacturer Contact

Austin , TX 78717

(512) 432 -9611

Device Event Key 331904

MDR Report Key 342566

Event Key 322596

Report Number 2935620-2001-01162

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

4/26/2004

Event Type Injury

Type of Report Initial

Report Date 06/18/2001

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 07/18/2001

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-049

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 06/18/2001

Was Device Evaluated By Manufacturer? No

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization Required Intervention

Adverse Event or Product Problem Description

Report Date: 06/18/2001 MDR Text Key: 1186749 Patient Sequence Number: 1

It was reported: metasul insert dislocation after 4 years implantation, significant metalosis. Severe pain.

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

# U.S. Food and Drug Administration

### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 11

New Search

Back To Search Results

## Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X53 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X53 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-053

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

SULZER ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Randy Jasek, Supervisor

9900 Spectrum Drive

Manufacturer Contact Austin, TX 78717

(512) 432 -9611

Device Event Key 341161

MDR Report Key 351880

**Event Key** 331492

Report Number 2935620-2001-01531

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Remedial Action Replace

Event Type No Answer Provided

Type of Report Initial

Report Date 08/23/2001

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 09/18/2001

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-053

Device LOT Number 1348667-B

Was Device Available For Evaluation? No

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 08/23/2001

Was Device Evaluated By Manufacturer? Device Not Returned To Manufacturer

Date Device Manufactured 08/01/2000

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization Required Intervention

Adverse Event or Product Problem Description
Report Date: 08/23/2001 MDR Text Key: 1221387 Patient Sequence Number: 1
It was reported: apr metasul insert disengaged [nine] months after surgery.

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

# FDA U.S. Food and Drug Administration & Department of Health and Human Services

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards

12

New Search

Back To Search Results

#### Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X53 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X53 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-053

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

SULZER ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Randy Jasek, Supervisor

9900 Spectrum Drive

Manufacturer Contact Austin, TX 78717

(512) 432 -9611

Device Event Key 355407

MDR Report Key 366284

**Event Key** 345365

Report Number 2935620-2001-01765

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

181

4/26/2004

Remedial Action Replace

Event Type Injury

Type of Report Initial

Report Date 11/26/2001

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 12/18/2001

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-053

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Date Returned to Manufacturer 11/26/2001

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 11/06/2001

Was Device Evaluated By Manufacturer? No

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization Required Intervention

Adverse Event or Product Problem Description

Report Date: 11/26/2001 MDR Text Key: 1272251 Patient Sequence Number: 1

It was reported: disassociation of liner from cup.

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page



FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 13

New Search

Back To Search Results

#### Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X49 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X49 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-049

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

SULZER ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Randy Jasek, Supervisor

9900 Spectrum Drive Manufacturer Contact

Austin , TX 78717

(512) 432 -9611

Device Event Key 359226

MDR Report Key 370104

**Event Key** 349096

Report Number 2935620-2001-01806

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Event Type Injury

Type of Report Initial

Report Date 12/12/2001

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 01/07/2002

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-049

Device LOT Number 1302992

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Date Returned to Manufacturer 12/12/2001

is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 12/12/2001

Was Device Evaluated By Manufacturer? No

Date Device Manufactured 06/01/1997

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization Required Intervention

Adverse Event or Product Problem Description

Report Date: 12/12/2001 MDR Text Key: 1285891 Patient Sequence Number: 1

It was reported: the anti-rotation pin was dislocated from the insert.

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

Center for Devices and Radiological Health / CDRH

4/26/2004

## FDA U.S. Food and Drug Administration & Department of Heelth and Human Services

#### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards

14

New Search

Back To Search Results

#### Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name FMP ACETABULAR SYSTEM

Type of Device ACETABULAR LINER

Baseline Brand Name FMP ACETABULAR SYSTEM

Baseline Generic Name ACETABULAR LINER

Baseline Catalogue Number 499-28-009

Baseline Device Family METAL/METAL LINER

Baseline Device 510(K) Number K003250

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 01/01/2001

ENCORE ORTHOPEDICS, INC.

Manufacturer (Section F) 9800 Metric Blvd.

Austin TX 78758

ENCORE ORTHOPEDICS, INC.

Manufacturer (Section D) 9800 Metric Blvd.

Austin TX 78758

9800 Metric Blvd.

Manufacturer Contact Austin, TX 78758

(800) 456 -8696

Device Event Key 359915

MDR Report Key 370788

Event Key 349781

Report Number 1644408-2002-00001

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Health Professional

279

4/26/2004

Event Type Other

Type of Report Initial

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 01/09/2002

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 499-28-009

Device LOT Number 770291

Was Device Available For Evaluation? Yes

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 12/17/2001

Was Device Evaluated By Manufacturer? Yes

Date Device Manufactured 10/01/2001

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Other

**Adverse Event or Product Problem Description** 

Report Date: 01/04/2001 MDR Text Key: 1288367 Patient Sequence Number: 1

Total hip replacement surgery extended due to failure to properly assemble acetabular liner into shell.

**Additional Manufacturer Narrative** 

Report Date: 01/09/2002 MDR Text Key: 1288368

Evaluation determined that event was caused by swelling of the polyethylene due to improper handling prior to surgery.

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

# U.S. Food and Drug Administration

## CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 15

New Search

Back To Search Results

## Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X51 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X51 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-051

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

SULZER ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Randy Jasek, Supervisor 9900 Spectrum Drive

Manufacturer Contact Austin , TX 78717

(512) 432 -9611

Device Event Key 364616

MDR Report Key 375510

**Event Key** 354337

Report Number 2935620-2001-01841

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

4/26/2004

Remedial Action Replace

Event Type Injury

Type of Report Initial

Report Date 01/28/2002

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 02/06/2002

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-051

Device LOT Number 1245308

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Date Returned to Manufacturer 01/28/2002

is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 12/05/2001

Was Device Evaluated By Manufacturer? No

Date Device Manufactured 06/01/1996

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization

Adverse Event or Product Problem Description

Report Date: 01/28/2002 MDR Text Key: 1305827 Patient Sequence Number: 1

It was reported dr. Reported he had a patient with pain. Disassociation of liner from shell.

**Patient TREATMENT DATA** 

Date Received: 02/06/2002 Patient Sequence Number: 1

Treatment Date Treatment

1 4311-00-051 APR II SCWLS SHELL 51MM (LOT#1295892)

01/01/2001

2,(2001)

3442 ALLOPRO BATORY FEM HD 28MM +4 (2001),01/01/2001

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 16

New Search

Back To Search Results

# Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X59 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X59 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-059

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

SULZER ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Allyson Hein, Qa Supervisor

9900 Spectrum Drive Manufacturer Contact

Austin, TX 78717

(512) 432 -9285

Device Event Key 367942

MDR Report Key 378913

Event Key 357627

Report Number 2935620-2002-00026

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Health Professional

Remedial Action Replace

Event Type Injury

Type of Report Initial

Report Date 01/28/2002

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 02/27/2002

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-059

Device LOT Number 1230120

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Date Returned to Manufacturer 01/28/2002

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 01/28/2002

Was Device Evaluated By Manufacturer? No

Date Device Manufactured 03/01/1996

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization Required Intervention

Adverse Event or Product Problem Description

Report Date: 01/28/2002 MDR Text Key: 1317161 Patient Sequence Number: 1

It was reported: disassociation of an apr metasul acetabular insert size 59mm. Patient fell in 2001.

#### Patient TREATMENT DATA

Date Received: 02/27/2002 Patient Sequence Number: 1

#		Treatment	Treatment Date
1	4311-00-059 APR II SHELL 59MM (LOT# 1298084)(2001)		01/01/2001

2 3461 AP BALL HD MED 28/14 (LOT# B023602)(2001). 01/01/2001

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer EDA Home Page | Search FDA Site | EDA A-Z Index | Contact EDA | HHS Home Page

### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 17

New Search

Back To Search Results

# Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X51 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X51 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-051

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

SULZER ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Allyson Hein, Supervisor

9900 Spectrum Drive Manufacturer Contact

Austin, TX 78717

(512) 432 -9285

Device Event Key 371588

MDR Report Key 382514

**Event Key** 361129

Report Number 2935620-2002-00069

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Health Professional

Remedial Action Replace

Event Type Injury

Type of Report Initial

Report Date 02/27/2002

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 03/19/2002

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-051

Device LOT Number, 1181946

Was Device Available For Evaluation? No

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Was Device Evaluated By Manufacturer? Device Not Returned To Manufacturer

Date Device Manufactured 09/01/1995

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization Required Intervention

Adverse Event or Product Problem Description

Report Date: 02/27/2002 MDR Text Key: 1329206 Patient Sequence Number: 1

It was reported: clinical advised on an apr metasul acetabular liner disassociation, patient was revised in 2002.

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

#### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 18

New Search

Back To Search Results

# Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X57 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28 X 57 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-057

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

SULZER ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Allyson Hein, Supervisor

9900 Spectrum Drive Manufacturer Contact

Austin , TX 78717

(512) 432 -9285

Device Event Key 376896

MDR Report Key 387867

**Event Key** 366275

Report Number 2935620-2002-00105

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Company Representative

Remedial Action Replace

Event Type Injury

Type of Report Initial

Report Date 03/15/2002

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 04/12/2002

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-057

Device LOT Number 1181949

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Date Returned to Manufacturer 03/28/2002

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 03/15/2002

Was Device Evaluated By Manufacturer? No

Date Device Manufactured 09/01/1995

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization Required Intervention

Adverse Event or Product Problem Description

Report Date: 03/15/2002 MDR Text Key: 1346855 Patient Sequence Number: 1

It was reported: disassocation of an apr metasul insert.

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search EDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

Center for Devices and Radiological Health / CDRH

#### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 19

New Search

Back To Search Results

### Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X51 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X51 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-051

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

SULZER ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Allyson Hein, Supervisor

9900 Spectrum Drive Manufacturer Contact Austin , TX 78717

(512) 432 -9285

Device Event Key 376891

MDR Report Key 387862

**Event Key** 366270

Report Number 2935620-2002-00106

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Company Representative

Remedial Action Replace

Event Type Injury

Type of Report Initial

Report Date 03/15/2002

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 04/12/2002

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-051

Device LOT Number 1245308

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Date Returned to Manufacturer 03/28/2002

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 03/15/2002

Was Device Evaluated By Manufacturer? No

Date Device Manufactured 06/01/1996

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization Required Intervention

Adverse Event or Product Problem Description

Report Date: 03/15/2002 MDR Text Key: 1346839 Patient Sequence Number: 1

It was reported: pain and breakage of implant.

Patient TREATMENT DATA

Date Received: 04/12/2002 Patient Sequence Number: 1

Treatment

Treatment Date

1 3441 ALLPRO BATTROY HD 28MM NEU (LOT#96583972)2002

01/01/2002

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page



#### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 20

**New Search** 

Back To Search Results

### Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X49 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X49 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-049

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

SULZER ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Allyson Hein, Supervisor

9900 Spectrum Drive

Manufacturer Contact Austin, TX 78717

(512) 432 -9285

Device Event Key 380445

MDR Report Key 391396

**Event Key** 369693

Report Number 2935620-2002-00141

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Company Representative

Remedial Action Replace

Event Type Injury

Type of Report Initial

Report Date 04/15/2002

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 05/01/2002

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-049

Device LOT Number 1322215

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Date Returned to Manufacturer 05/01/2002

is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 04/15/2002

Was Device Evaluated By Manufacturer? No

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization Required Intervention

**Adverse Event or Product Problem Description** 

Report Date: 04/15/2002 MDR Text Key: 1359227 Patient Sequence Number: 1

It was reported: disassociation of insert 2 years 9 months after initial surgery.

Patient TREATMENT DATA

Date Received: 05/01/2002 Patient Sequence Number: 1

**Treatment Date** Treatment

01/01/2002

1 3460 METASUL BALL HD 28MM SZ S (LOT#B436274)(2002)

Database contains data received through March 30, 2004

CDRH Home Page | CDRH.A.-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search EDA Site | FDA A.-Z Index | Contact FDA | HHS Home Page

#### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 21

New Search

Back To Search Results

### Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name COCR HD 12/14 NECK NEU 28MM METASUL

Type of Device HIP PROSTHESIS

Baseline Brand Name COCR HD 12/14 NECK NEU 28MM METASUL

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 7340-28-000

Baseline Device Family NA

Baseline Device 510(K) Number K974728

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 08/03/1999

SULZER ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

SULZER ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Allyson Hein, Supervisor

9900 Spectrum Drive Manufacturer Contact

Austin , TX 78717

(512) 432 -9285

Device Event Key 390758

MDR Report Key 401722

Event Key 379644

Report Number 2935620-2002-00178

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Company Representative

Remedial Action Replace

Event Type Injury

Type of Report Initial

Report Date 06/04/2002

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 06/26/2002

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 7340-28-000

Device LOT Number 1426357

Was Device Available For Evaluation? No

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 06/04/2002

Was Device Evaluated By Manufacturer? Device Not Returned To Manufacturer

Date Device Manufactured 07/01/2000

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization Required Intervention

Adverse Event or Product Problem Description

Report Date: 06/04/2002 MDR Text Key: 1395451 Patient Sequence Number: 1

It was reported: patient was revised due to dislocation of the hip.

Database contains data received through March 30, 2004

CDRH Home Page | CDRH.A.Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A.Z Index | Contact FDA | HHS Home Page



CENTER FOR DEVICES AND RADIOLOGICAL HEALTH FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 22

New Search

Back To Search Results

### Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name M2A MODULAR HEAD

Type of Device PROSTHESIS, HIP, COMP.

Baseline Brand Name M2A MODULAR HEAD

Baseline Generic Name PROSTHESIS, HIP, COMP

Baseline Catalogue Number 11-173660

Baseline Device Family M2A MODULAR HEAD

Baseline Device 510(K) Number K011110

Baseline Shelf Life Information Yes

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) 12

Date First Marketed 07/02/2001

BIOMET, INC.

Manufacturer (Section F) P.O. Box 587

Warsaw IN 46581 0587

BIOMET, INC.

Manufacturer (Section D) P.O. Box 587

Warsaw IN 46581 0587

Beth Albert, Correct Action As P.O. Box 587

Manufacturer Contact

Warsaw, IN 46581-0587

(574) 267 -6639

Device Event Key 417045

MDR Report Key 428854

**Event Key** 404930

Report Number 1825034-2002-00129

Device Sequence Number 1

Product Code KWA

4/26/2004

Report Source Manufacturer

Source Type User facility, Company Representative

Event Type Injury

Type of Report Initial

Report Date 11/08/2002

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 11/20/2002

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device EXPIRATION Date 07/01/2012

Device Catalogue Number 11-173660

Device LOT Number 698740

Was Device Available For Evaluation? No

Date Manufacturer Received 11/08/2002

Was Device Evaluated By Manufacturer? Device Not Returned To Manufacturer

Date Device Manufactured 07/01/2002

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization Required Intervention

#### Adverse Event or Product Problem Description

Report Date: 11/08/2002 MDR Text Key: 1488471 Patient Sequence Number: 1

Report rec'd from hospital stating that during the time period from event date to 11/2002, eleven total joint and/or revision procedures were performed, with five of these cases resulting in post operative infection. The procedures performed are comprised of various components being implanted. The co has no info that these incidents of infections are in any way related to the biomet device. The co is continuing to gather further info, but to date there is no common thread involving the device or its components that would lead biomet to conclude that the device caused any infection. Investigation continues at the hospital for non-device causes. Left "tha" performed in 2002, resulting in postoperative infection.

#### Additional Manufacturer Narrative

Report Date: 11/08/2002 MDR Text Key: 1488474

Biomet does not believe that the company's devices are implicated in these infections because similar reports of infection have not been rec'd from other user facilities for any of the listed mfg lots. Notwithstanding the fact that the company does not believe these events are related, the company is advising fda of these incidents under part 803 so that fda's files are complete. It is the co understanding that the institution has filed a medwatch form for these infections and biomet wants fda to know the status of its investigation.

Database contains data received through March 30, 2004

# U.S. Food and Drug Administration CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 23

New Search

Back To Search Results

# Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name COCR HD 12/14 +4MM NECK 28MM METASUL

Type of Device HIP PROSTHESIS

Baseline Brand Name COCR HD 12/14 +4MM NECK 28MM METASUL

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 7340-28-400

Baseline Device Family NA

Baseline Device 510(K) Number K974728

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 08/03/1999

CENTERPULSE ORTHOPEDICS, INC

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

CENTERPULSE ORTHOPEDICS, INC

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Ron Yarbrough, Supervisor

9900 Spectrum Dr

Manufacturer Contact Austin, TX 78717

(512) 432 -9437

Device Event Key 421879

MDR Report Key 432921

Event Key 409618

Report Number 2935620-2002-00393

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Company Representative

Remedial Action Replace

Event Type Injury

Type of Report Initial

Report Date 12/10/2002

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 12/17/2002

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 7340-28-400

Device LOT Number 1409774

Was Device Available For Evaluation? No

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 12/10/2002

Was Device Evaluated By Manufacturer? Device Not Returned To Manufacturer

Date Device Manufactured 10/01/1999

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization Required Intervention

**Adverse Event or Product Problem Description** 

Report Date: 12/10/2002 MDR Text Key: 1502746 Patient Sequence Number: 1

It was reported: 2nd revision was required because of the instability of the hip, and was operated in 2001. This is apparent, given the fact that they converted to a hooded insert, as well as the fact that a +8 cocr head was used. Both of these changes will help to increase the stability of the hip. 1st rev. Of this patient is of a recall shell. P01-2916 (rm database).

Database contains data received through March 30, 2004

CDRH Home Page | CDRH.A.-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A.-Z Index | Contact FDA | HHS Home Page

Center for Devices and Radiological Health / CDRH

4/26/2004

# CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 24

New Search

Back To Search Results

# Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name M2A MODULAR HEAD

Type of Device PROSTHESIS, HIP, COMP.

Baseline Brand Name M2A MODULAR HEAD

Baseline Generic Name PROSTHESIS, HIP, COMP

Baseline Catalogue Number 11-173661

Baseline Device Family M2A MODULAR HEAD

Baseline Device 510(K) Number K011110

Baseline Shelf Life Information Yes

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) 12

Date First Marketed 07/02/2001

BIOMET, INC.

Manufacturer (Section F) P.O. Box 587

Warsaw IN 46581 0587

BIOMET, INC.

Manufacturer (Section D) P.O. Box 587

Warsaw IN 46581 0587

Beth Albert, Correct Action Ast

P.O. Box 587

Manufacturer Contact Warsaw, IN 46581-0587

(574) 267 -6639

Device Event Key 425634

MDR Report Key 436691

**Event Key** 413272

Report Number 1825034-2003-00001

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Company Representative

Event Type Injury

Type of Report Initial

Report Date 01/07/2003,12/09/2002

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 01/07/2003

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device EXPIRATION Date 04/01/2012

Device Catalogue Number 11-173661

Device LOT Number 321210

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Date Returned to Manufacturer 12/09/2002

Is The Reporter A Health Professional? Yes

Was the Report Sent to FDA? No

Distributor Facility Aware Date 12/10/2002

Device Age 8 mo

Event Location Hospital

Date Manufacturer Received 12/09/2002

Was Device Evaluated By Manufacturer? Yes

Date Device Manufactured 04/01/2002

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization Required Intervention

Adverse Event or Product Problem Description

Report Date: 01/07/2003 MDR Text Key: 1515315 Patient Sequence Number: 1

Following total hip arthroplasty performed in 2002, pt continued to experience hip pain. Pt underwent add'l surgery 1 mo later, where osteophytes were removed and modular head component was exchanged.

Additional Manufacturer Narrative

Report Date: 01/07/2003 MDR Text Key: 1515318

There are warnings in the package insert that state that this type of event can occur. This type of event is not occurring at a rate above expected frequency. No remedial action will be taken. No further complications have been reported. Current info is insufficient to permit a valid conclusion as to the cause of the event.

207

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 25

New Search

Back To Search Results

# Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ULTAMET MOM INSRT 540DX36ID

Type of Device TOTAL HIP PROSTHESIS

Baseline Brand Name ULTAMET MOM INSRT 540DX36ID

Baseline Generic Name PINNACLE METAL ON METAL ACETABULAR CUP

Baseline Catalogue Number 121887354

Baseline Device Family ULTAMET MOM INSERT

Baseline Device 510(K) Number K003523

Baseline Shelf Life Information No

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/13/2000

DEPUY INTERNATIONAL,LTD.

St. Anthony's Rd Manufacturer (Section F) Beeston, Leeds

UNITED KINGDOM LS11 8DT

DEPUY INTERNATIONAL, LTD.

St. Anthony's Rd Manufacturer (Section D) Beeston, Leeds

UNITED KINGDOM LS11 8DT

DEPUY ORTHOPAEDICS, INC.

Manufacturer (Section G) 700 Orthopaedic Drive

Warsaw IN 46581 0988

Hans Kusserow, Mgr.

700 Orthopaedic Drive Manufacturer Contact Warsaw , IN 46581-0988

(574) 372 -7416

Device Event Key 426244

MDR Report Key 437308

**Event Key** 413857

Report Number 1818910-2003-00019

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Distributor

Remedial Action Other

Event Type Malfunction

Type of Report Initial

Report Date 01/10/2003

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 01/10/2003

Is This An Adverse Event Report? No

Is This A Product Problem Report? Yes

Device Operator Health Professional

Device Catalogue Number 121887354

Device LOT Number YMD59

Was Device Available For Evaluation? Yes

Date Returned to Manufacturer 12/13/2002

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Distributor Facility Aware Date 12/13/2002

Device Age unknown

Event Location Hospital

Date Manufacturer Received 12/11/2002

Was Device Evaluated By Manufacturer? Yes

Date Device Manufactured 07/01/2002

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Other

Adverse Event or Product Problem Description

Report Date: 01/10/2003 MDR Text Key: 1517381 Patient Sequence Number: 1

The package is split - inner plastic portion.

Additional Manufacturer Narrative

Report Date: 01/10/2003 MDR Text Key: 1517384

This complaint is still under investigation. Depuy will notify the fda of the results of this investigation once it has been completed.

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

#### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 26

New Search

Back To Search Results

# Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ULTAMET MOM INSRT 520DX36ID

Type of Device TOTAL HIP PROSTHESIS

Baseline Brand Name ULTAMET MOM INSERT 520DX36ID

Baseline Generic Name PINNACLE METAL ON METAL ACETABULAR CUP

Baseline Catalogue Number 121887352

Baseline Device Family ULTAMET MOM INSERT

Baseline Device 510(K) Number K003523

Baseline Shelf Life Information No

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/13/2000

DEPUY INTERNATIONAL, LTD.

St. Anthony's Road Manufacturer (Section F)

Leeds

UNITED KINGDOM LS11 8DT

DEPUY INTERNATIONAL, LTD.

St. Anthony's Road Manufacturer (Section D)

Leeds

UNITED KINGDOM LS11 8DT

DEPUY ORTHOPAEDICS, INC.

Manufacturer (Section G) 700 Orthopaedic Drive

Warsaw IN 46581 0988

Hans Kusserow, Mgr.

700 Orthopaedic Drive Manufacturer Contact Warsaw , IN 46581-0988

(574) 372 -7416

Device Event Key 439175

MDR Report Key 450176

**Event Key** 426312

306

ลเฉ

Report Number 1818910-2003-00178

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Health Professional, Distributor

Remedial Action Other

Event Type Injury

Type of Report Initial

Report Date 03/27/2003

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 03/27/2003

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? Yes

Device Operator Health Professional

Device Catalogue Number 121887352

Device LOT Number 1070181

Was Device Available For Evaluation? Yes

Date Returned to Manufacturer 03/03/2003

Is The Reporter A Health Professional? Yes

Was the Report Sent to FDA? No

Distributor Facility Aware Date 02/26/2003

Device Age unknown

Event Location Hospital

Date Manufacturer Received 02/26/2003

Was Device Evaluated By Manufacturer? Yes

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Required Intervention

Adverse Event or Product Problem Description

Report Date: 03/27/2003 MDR Text Key: 1562106 Patient Sequence Number: 1

The metal insert was locked into the shell off-axis. Insert could not be removed from acetabular shell.

Additional Manufacturer Narrative

Report Date: 03/27/2003 MDR Text Key: 1562109

This complaint is still under investigation. Depuy will notify the fda of the results of this investigation once it has been completed.

CORH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

# U.S. Food and Drug Administration CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 27

New Search

Back To Search Results

### Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X49 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X49 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-049

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

CENTERPULSE ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

CENTERPULSE ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

Jim Gonzales, Supervisor

9900 Spectrum Drive Manufacturer Contact

Austin , TX 78717

(512) 432 -9688

Device Event Key 456610

MDR Report Key 467685

**Event Key** 443229

Report Number 2935620-2003-00135

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Company Representative

Remedial Action Replace

Event Type Injury

Type of Report Initial

Report Date 06/10/2003

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 06/25/2003

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-049

Device LOT Number 1230114

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Date Returned to Manufacturer 06/13/2003

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 06/10/2003

Was Device Evaluated By Manufacturer? No

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization

Adverse Event or Product Problem Description Report Date: 06/10/2003 MDR Text Key: 1622790 Patient Sequence Number: 1 Pt felt pain, due to dislocation. Pt was revised.

#### Patient TREATMENT DATA

Date Received: 06/25/2003 Patient Sequence Number: 1

# Treatment Treatment Date

1 4310-02-049 APR 12 SLOT SHELL SZ49MM (LOT#1268548)

01/01/2003

2 3462 AP BALL HEAD LONG 28/14 (LOT#96547339).

01/01/2003

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disciaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 28

New Search

Back To Search Results

# Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X51 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X51 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-051

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

CENTERPULSE ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

CENTERPULSE ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

CENTERPULSE ORTHOPEDICS, INC.

Manufacturer (Section G) 9900 Spectrum Drive

Austin TX 78717

Jim Gonzales, Supervisor

9900 Spectrum Drive Manufacturer Contact Austin , TX 78717

(512) 432 -9688

Device Event Key 460922

MDR Report Key 472071

**Event Key** 447449

Report Number 2935620-2003-00143

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Company Representative

Remedial Action Replace

Event Type Injury

Type of Report Initial

Report Date 06/19/2003

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 07/16/2003

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-051

Device LOT Number 1303680

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Date Returned to Manufacturer 06/20/2003

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 06/19/2003

Was Device Evaluated By Manufacturer? No

Date Device Manufactured 11/01/1997

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization

Adverse Event or Product Problem Description

Report Date: 06/19/2003 MDR Text Key: 1637251 Patient Sequence Number: 1

It was reported: dislocation of apr metasul insert.

Patient TREATMENT DATA

Date Received: 07/16/2003 Patient Sequence Number: 1

Treatment

1 3461 AP BALL HEAD MED 28/14 (LOT#B525466)(2002).

**Treatment Date** 

01/01/2002

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

#### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 29

New Search

Back To Search Results

### Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X49 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X49 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-049

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

CENTERPULSE ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

CENTERPULSE ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

CENTERPULSE ORTHOPEDICS, INC.

Manufacturer (Section G) 9900 Spectrum Dr

Austin TX 78717

Jim Gonzales, Supervisor

9900 Spectrum Drive Manufacturer Contact

Austin , TX 78717 (512) 432 -9688

Device Event Key 461150

MDR Report Key 472298

**Event Key** 447676

Report Number 2935620-2003-00142

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Company Representative

Remedial Action Replace

Event Type Injury

Type of Report Initial

Report Date 06/25/2003

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 07/17/2003

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-049

Device LOT Number 1322216

Was Device Available For Evaluation? No

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 06/25/2003

Was Device Evaluated By Manufacturer? Device Not Returned To Manufacturer

Date Device Manufactured 05/01/1998

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization

**Adverse Event or Product Problem Description** 

Report Date: 06/25/2003 MDR Text Key: 1638072 Patient Sequence Number: 1

Pt had a dislocation and was revised.

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 30

New Search

Back To Search Results

# Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name SZ 28X55MM STD METASUL INS I-O

Type of Device HIP PROSTHESIS

Baseline Brand Name SZ 28 X 55MM STD METASUL INS I-O

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4372-28-055

Baseline Device Family NA

Baseline Device 510(K) Number K974728

ls Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 08/03/1999

CENTERPULSE ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

CENTERPULSE ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

CENTERPULSE ORTHOPEDICS, INC.

Manufacturer (Section G) 9900 Spectrum Dr

Austin TX 78717

Jim Gonzales, Supervisor

9900 Spectrum Drive Manufacturer Contact Austin , TX 78717

(512) 432 -9688

Device Event Key 478775

MDR Report Key 490033

Event Key 464565

Report Number 2935620-2003-00246

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Company Representative

Remedial Action Other

Event Type Other

Type of Report Initial

Report Date 10/14/2003

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 10/14/2003

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4372-28-055

Device LOT Number 1426338-A

Was Device Available For Evaluation? No

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 10/09/2003

Was Device Evaluated By Manufacturer? No

Date Device Manufactured 08/01/2002

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization

Adverse Event or Product Problem Description Report Date: 10/14/2003 MDR Text Key: 1696369 Patient Sequence Number: 1 Delay in surgery.

Database contains data received through March 30, 2004

CORH Home Page | CORH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

#### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 31

New Search

Back To Search Results

### Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X49 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X49 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-049

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

CENTERPULSE ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

CENTERPULSE ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

CENTERPULSE ORTHOPEDICS, INC.

Manufacturer (Section G) 9900 Spectrum Dr

Austin TX 78717

Jim Gonzales, Supervisor

9900 Spectrum Drive Manufacturer Contact Austin, TX 78717

(512) 432 -9688

Device Event Key 489867

MDR Report Key 501097

Event Key 475206

Report Number 2935620-2003-00279

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Company Representative

Remedial Action Replace

Event Type No Answer Provided

Type of Report Initial

Report Date 12/16/2003

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 12/16/2003

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-049

Device LOT Number 1230114

Was Device Available For Evaluation? No

is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 11/26/2003

Was Device Evaluated By Manufacturer? Device Not Returned To Manufacturer

Date Device Manufactured 03/01/1996

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization

Adverse Event or Product Problem Description Report Date: 12/16/2003 MDR Text Key: 1732417 Patient Sequence Number: 1

It was reported: pt was revised.

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

## U.S. Food and Drug Administration & Department of Health and Human Services

### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 32

**New Search** 

Back To Search Results

## Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X55 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X55 METASUL APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-055

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

CENTERPULSE ORTHOPEDICS, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

CENTERPULSE ORTHOPEDICS, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

CENTERPULSE ORTHOPEDICS, INC.

Manufacturer (Section G) 9900 Spectrum Drive

Austin TX 78717

Jim Gonzales, Supervisor

Manufacturer Contact 9900 Spectrum Dr Austin . TX 78717

(512) 432 -9688

Device Event Key 491229

MDR Report Key 502457

**Event Key** 476487

Report Number 2935620-2003-00282

Device Sequence Number 1

320

226

Product Code KWA

Report Source Manufacturer

Source Type Company Representative

Remedial Action Replace

Event Type Injury

Type of Report Initial

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 12/23/2003

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-055

Was Device Available For Evaluation? Yes

Date Returned to Manufacturer 12/17/2003

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 12/17/2003

Was Device Evaluated By Manufacturer? No

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization

Adverse Event or Product Problem Description

Report Date: 10/23/2003 MDR Text Key: 1737004 Patient Sequence Number: 1

It was reported: patient was revised in 2003.

Database contains data received through March 30, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

# U.S. Food and Drug Administration

### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title 21 | Advisory Committees | Assembler | NHRIC | Guidance | Standards 33

**New Search** 

Back To Search Results

## Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ULTAMET MOM INSRT 620DX36ID

Type of Device TOTAL HIP REPLACEMENT

Baseline Brand Name ULTAMET MOM INSERT 620DX36ID

Baseline Generic Name PINNACLE METAL ON METAL ACETABULAR CUP

Baseline Catalogue Number 121887362

Baseline Device Family ULTAMET MOM INSERT

Baseline Device 510(K) Number K003523

Baseline Shelf Life Information No

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 12/13/2000

DEPUY ORTHOPAEDICS, INC.

700 Orthopaedic Dr.

Manufacturer (Section F) P.O. Box 988

Warsaw IN 46581 0988

DEPUY ORTHOPAEDICS, INC.

700 Orthopaedic Dr.

Manufacturer (Section D) P.O. Box 988

Warsaw IN 46581 0988

DEPUY ORTHOPAEDICS, INC.

Manufacturer (Section G) 700 Prthopaedic Drive

Warsaw IN 46581 0988

Hans Kusserow, Mgr.

700 Orthopaedic Drive Manufacturer Contact Warsaw , IN 46581-0988

(574) 372 -7416

Device Event Key 495817

MDR Report Key 506895

**Event Key** 480781

322

http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfMAUDE/Detail.CFM?MDRFOI\_\_ID=506895

Report Number 1818910-2004-00044

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Health Professional, Distributor

Remedial Action Other

Event Type Injury

Type of Report Initial

Report Date 12/19/2003

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 01/16/2004

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? Yes

Device Operator Health Professional

Device Catalogue Number 121887362

Device LOT Number YHT84

Was Device Available For Evaluation? Yes

Is The Reporter A Health Professional? Yes

Was the Report Sent to FDA? No

**Event Location** Hospital

Date Manufacturer Received 12/19/2003

Was Device Evaluated By Manufacturer? Yes

Date Device Manufactured 03/01/2002

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Required Intervention

Adverse Event or Product Problem Description

Report Date: 12/19/2003 MDR Text Key: 1751674 Patient Sequence Number: 1

Metal insert did not seat properly. Surgery was extended by 30 minutes due to the product problem.

Additional Manufacturer Narrative

Report Date: 12/19/2003 MDR Text Key: 1751677

This complaint is still under investigation. Depuy will notify the fda of the results of this investigation

once it has been completed.

Database contains data received through March 30, 2004

34

## U.S. Food and Drug Administration & Department of Health and Human Services

#### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title | Advisory | Assembler | NHRIC | Guidance | Standards 21 | Committees

New Search

Back To Search Results

### Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name SROM INS M28 10DEG 28MM C

Type of Device TOTAL HIP PROSTHESIS

Baseline Brand Name SROM INS M28 10DEG 28MM C

Baseline Generic Name SROM INSERT

Baseline Catalogue Number 552710

Baseline Device Family SROM INSERT 10 DEGREE

Baseline Device 510(K) Number K924492

Baseline Shelf Life Information No

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 11/24/1992

DEPUY ORTHOPAEDICS, INC.

Manufacturer (Section F) 700 Orthopaedic Dr.

Warsaw IN 46581 0988

DEPUY ORTHOPAEDICS, INC.

Manufacturer (Section D) 700 Orthopaedic Dr.

Warsaw IN 46581 0988

DEPUY ORTHOPAEDICS, INC.

Manufacturer (Section G) 700 Orthopaedic Dr

Warsaw IN 46581 0986

Hans Kusserow, Mgr. 700 Orthopaedic Drive

Manufacturer Contact Warsaw, IN 46581-0988

(574) 372 -7416

Device Event Key 494588

231

MDR Report Key 505719

**Event Key** 479649

Report Number 1818910-2004-00050

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Health Professional, Distributor

Remedial Action Other

Event Type Injury

Type of Report Initial

Report Date 12/10/2003

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 01/09/2004

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? Yes

**Device Operator** Health Professional

**Device Catalogue Number** 552710

Was Device Available For Evaluation? Yes

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Distributor Facility Aware Date 12/10/2003

Device Age unknown

**Event Location** Hospital

Date Manufacturer Received 12/10/2003

Was Device Evaluated By Manufacturer? Yes

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Required Intervention

Adverse Event or Product Problem Description

Report Date: 12/10/2003 MDR Text Key: 1747829 Patient Sequence Number: 1

Ring broken, worn liner. Revised liner only to non-constrained.

**Additional Manufacturer Narrative** 

Report Date: 12/10/2003 MDR Text Key: 1747832

The investigation could not firmly establish a root cause for the event. Based on the investigation, this is the first reported complaint of wear or breakage against this

product code, and as such, the need for corrective action has not been indicated. Should additional info related to this event be rec'd, the investigation will re-open. Depuy considers the investigation closed at this time.

Database contains data received through June 29, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page



#### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title | Advisory | Assembler | NHRIC | Guidance | Standards

35

**New Search** 

Back To Search Results

### Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name S-ROM INSERT L32, 10DEG., 32MM C

Type of Device TOTAL HIP REPLACEMENT

Baseline Brand Name S-ROM INSERT L32, 10 DEG., 32MM C

Baseline Generic Name SROM INSERT

Committees

Baseline Catalogue Number 521632

Baseline Device Family SROM INSERT

Baseline Device 510(K) Number K924492

Is Baseline PMA Number Provided?

Baseline Preamendment? No

Manufacturer (Section F)

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 11/24/1992

DEPUY-RAYNHAM A DIV OF DEPUY

ORTHOPAEDICS, INC.

anulacturer (Section F) 325 Paramount Dr

Raynham MA 02767 0350

DEPUY-RAYNHAM A DIV OF DEPUY

Manufacturer (Section D) ORTHOPAEDICS, INC.

325 Paramount Dr Raynham MA 02767 0350

DEPUY ORTHOPAEDICS, INC.

Manufacturer (Section G) 325 Paramount Dr

Raynham MA 02767 0350

Hans Kusserow, Mgr.
700 Orthopaedic Drive

Warsaw , IN 46581-0988

(574) 372 -7416

234

**Device Event Key** 501622

MDR Report Key 512613

**Event Key** 486353

Report Number 1818910-2004-00147

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Health Professional, Other

Remedial Action Other

Event Type Injury

Type of Report Initial

Report Date 02/10/2004

1 Device Was Involved in the

1 Patient Was Involved in the

Event

Date FDA Received 02/24/2004

Is This An Adverse Event Yes

Report?

Is This A Product Problem Yes

Report?

**Device Operator** Health Professional

**Device Catalogue Number** 521632

Device LOT Number SSC100739

Was Device Available For No

**Evaluation?** 

Is The Reporter A Health

Professional?

Was the Report Sent to FDA? No

Distributor Facility Aware Date 02/10/2004

Device Age 7 yr

**Event Location** Hospital

Date Manufacturer Received 02/10/2004

Was Device Evaluated By Manufacturer? Device Not Returned To Manufacturer

Date Device Manufactured 04/01/1994

Is The Device Single Use? Yes

Type of Device Usage Initial

#### Patient Outcome Required Intervention

Adverse Event or Product Problem Description

Report Date: 02/10/2004 MDR Text Key: 1770308 Patient Sequence Number: 1

Attorney reports this pt was revised due to premature polyethylene failure.

**Additional Manufacturer Narrative** 

Report Date: 02/10/2004 MDR Text Key: 1770311

This complaint is still under investigation. Depuy will notify the fda of the results of

this investigation once it has been completed.

Database contains data received through June 29, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

Committees

## FDA U.S. Food and Drug Administration & Department of Health and Human Services

#### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA CFR Title | Advisory | Assembler | NHRIC | Guidance | Standards 36

New Search

Back To Search Results

## Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ALPHA MTSL 28 STD JJ/28

Type of Device HIP PROSTHESIS

Baseline Brand Name ALPHA MTSL 28 STD JJ/28

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 01.00010.410

Baseline Device Family NA

Baseline Device 510(K) Number K003758

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

510(K) Exempt? No

Shelf Life(Months) NA

Date First Marketed 03/07/2001

CENTERPULSE ORTHOPEDICS

LTD.

Manufacturer (Section F) Po Box Ch-8404

Winterthur

SWITZERLAND

CENTERPULSE ORTHOPEDICS

LTD.

Manufacturer (Section D) Po Box Ch-8404

Winterthur

**SWITZERLAND** 

Ron Yarbrough, Supervisor

Manufacturer Contact 9900 Spectrum Drive

Austin , TX 78717

(512) 432 -9437

Device Event Key 434498

237

MDR Report Key 445528

**Event Key** 421802

Report Number 9613350-2003-00003

Device Sequence Number 1

Product Code KWA

Report Source Manufacturer

Source Type Company Representative

Remedial Action Replace

Event Type Injury

Type of Report Initial

Report Date 02/13/2003

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 02/28/2003

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 01.00010.410

Device LOT Number B661398

Was Device Available For Evaluation? No

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 02/13/2003

Was Device Evaluated By Device Not Returned To Manufacturer Manufacturer?

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization Required Intervention

Adverse Event or Product Problem Description

Report Date: 02/13/2003 MDR Text Key: 1545615 Patient Sequence Number: 1

Revision.

Additional Manufacturer Narrative

Report Date: 02/13/2003 MDR Text Key: 1545618

This investigation is considered closed. If any pertinent information becomes available this investigation will be reevaluated at that time.

Database contains data received through June 29, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page



#### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



 510(k)
 Registration
 Listing
 Adverse Events
 PMA
 Classification
 CLIA

 CFR Title
 | Advisory
 | Assembler
 | NHRIC
 | Guidance
 | Standards

 21
 Committees

37

New Search

Back To Search Results

#### Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name S-ROM INSERT M26, 15DEG, 26MM

Type of Device TOTAL HIP REPLACEMENT

Baseline Brand Name S-ROM INSERT M28, 15DEG, 26MM

Baseline Generic Name S-ROM INSERT

Baseline Catalogue Number 551946

Baseline Device Family S-ROM INSERT

Baseline Device 510(K) Number K924492

Baseline Shelf Life Information No

Is Baseline PMA Number Provided? No

Baseline Preamendment? No

Transitional? No

**510(K) Exempt?** No

Shelf Life(Months) NA

Date First Marketed 11/24/1992

DEPUY ORTHOPAEDICS, INC.

**Manufacturer** (Section F) 700 Orthopaedic Dr.

Warsaw IN 46581 0988

DEPUY ORTHOPAEDICS, INC.

Manufacturer (Section D) 700 Orthopaedic Dr.

Warsaw IN 46581 0988

DEPUY ORTHOPAEDICS, INC.

Manufacturer (Section G) 700 Orthopaedics Drive

Warsaw IN 46581 0988

Hans Kusserow, Mgr.

Manufacturer Contact 700 Orthopaedic Drive Warsaw, IN 46581-0988

(574) 372 -7416

Device Event Key 515787

240

MDR Report Key 526619

**Event Key** 499874

Report Number 1818910-2004-00385

**Device Sequence Number 1** 

Product Code KWA

Report Source Manufacturer

Source Type Distributor

Remedial Action Other

**Event Type** Malfunction

Type of Report Initial

**Report Date** 05/19/2004

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 05/25/2004

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

**Device Operator** Health Professional

**Device Catalogue Number** 551946

Device LOT Number SC104521

Was Device Available For Evaluation? No

Is The Reporter A Health Professional? No Answer Provided

Was the Report Sent to FDA? No

**Distributor Facility Aware Date** 05/19/2004

Device Age unknown

**Event Location** Hospital

Date Manufacturer Received 05/19/2004

Was Device Evaluated By Manufacturer? Device Not Returned To Manufacturer

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Required Intervention

**Adverse Event or Product Problem Description** 

Report Date: 05/19/2004 MDR Text Key: 1818260 Patient Sequence Number: 1 Surgeon claims polyethylene wear started to appear on x-rays 3 years after tha.

**Additional Manufacturer Narrative** 

Report Date: 05/19/2004 MDR Text Key: 1818263

This complaint is still under investigation. Depuy will notify the fda of the results of

this investigation once it has been completed. This complaint originated in japan. Therefore, we only have the hospital name.

Database contains data received through June 29, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

Committees

## U.S. Food and Drug Administration & Department of Health and Human Services

#### CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

FDA Home Page | CDRH Home Page | Search | CDRH A-Z Index | Contact CDRH



 510(k) | Registration | Listing | Adverse Events | PMA | Classification | CLIA

 CFR Title | Advisory | Assembler | NHRIC | Guidance | Standards

38

New Search

Back To Search Results

### Manufacturer and User Facility Device Experience (MAUDE) Database

Brand Name ACET INS SZ 28X49 METASUL-APR

Type of Device HIP PROSTHESIS

Baseline Brand Name ACET INS SZ 28X49 METASUL-APR

Baseline Generic Name HIP PROSTHESIS

Baseline Catalogue Number 4340-28-049

Baseline Device Family NA

Baseline Device 510(K) Number K993569

Is Baseline PMA Number Provided? No

**Baseline Preamendment?** No

Transitional? No

**510(K) Exempt?** No

Shelf Life(Months) NA

Date First Marketed 12/01/1999

ZIMMER AUSTIN, INC.

Manufacturer (Section F) 9900 Spectrum Dr.

Austin TX 78717

ZIMMER AUSTIN, INC.

Manufacturer (Section D) 9900 Spectrum Dr.

Austin TX 78717

ZIMMER AUSTIN, INC.

Manufacturer (Section G) 9900 Spectrum Dr

Austin TX 78717

Jim Gonzales, Supervisor

Manufacturer Contact 9900 Spectrum Dr

Austin, TX 78717 (512) 432 -9688

Device Event Key 520393

MDR Report Key 531144

**Event Key** 504236

Report Number 2935620-2004-00061

**Device Sequence Number 1** 

**Product Code KWA** 

Report Source Manufacturer

Source Type Company Representative

Remedial Action Replace

Event Type Injury

Type of Report Initial

**Report Date** 06/24/2004

1 Device Was Involved in the Event

1 Patient Was Involved in the Event

Date FDA Received 06/24/2004

Is This An Adverse Event Report? Yes

Is This A Product Problem Report? No

Device Operator Health Professional

Device Catalogue Number 4340-28-049

**Device LOT Number** 11/87758

Was Device Available For Evaluation? Device Not Returned To Manufacturer

Date Returned to Manufacturer 06/21/2004

Is The Reporter A Health Professional? No

Was the Report Sent to FDA? No

Date Manufacturer Received 06/21/2004

Was Device Evaluated By Manufacturer? No

Date Device Manufactured 09/01/1995

Is The Device Single Use? Yes

Type of Device Usage Initial

Patient Outcome Hospitalization

Adverse Event or Product Problem Description

Report Date: 06/24/2004 MDR Text Key: 1834094 Patient Sequence Number: 1

It was reported: pt was revised in 2004.

Patient TREATMENT DATA

Date Received: 06/24/2004 Patient Sequence Number: 1

# Treatment Treatment Date

1,4310-02-049 12 SLOT SHELL 49MM APR II(LOT#1202344),

2,(2004), 3,3461 AP BALL HEAD MED 28/14(LOT#95336363), 4,(2004).,

Database contains data received through June 29, 2004

CDRH Home Page | CDRH A-Z Index | Contact CDRH | Accessibility | Disclaimer FDA Home Page | Search FDA Site | FDA A-Z Index | Contact FDA | HHS Home Page

#### Class III Device Summary

#### Metal on Metal Articulating Surfaces

#### II. Literature Review of Risks Associated with Metal on Metal Articulating Surfaces

The following is a list of types of safety and effectiveness concerns/issues associated with the metal on metal total hip arthroplasty reported in published literature. Refer to following pages for the Bibliography of Abstracts, followed by copies of the abstract.

- 1. Radiological changes
- 2. Fretting corrosion
- 3. Aseptic loosening
- 4. Elevated serum levels of Cobalt and/or Chromium (significance unknown)
- 5 Pain
- 6. Dislocation
- 7. Revision
- 8. Septic loosening
- 9. Periarticular ossification
- 10. Elevated erythrocytes and urine metal ions
- 11. Early loosening
- 12. Infection
- 13. Calcar resorption
- 14. Artifact arising from metal hardware
- 15. Ectopic ossification
- 16. Progressive cement/bone interface radiolucencies
- 17. Cancer risk (increased risk not established)
- 18. Metal sensitivity
- 19. Metal wear debris

## METAL ON METAL TOTAL HIP ARTHROPLASTY BIBLIOGRAPHY OF ABSTRACTS

- Adami, G., Smarrelli, D., Martinelli, B., Acquavita, A., Reisenhofer, E.. "Cobalt Blood Levels After Total Hip Replacement (THR): A New Follow-Up Study In Trieste (Italy)". <u>Ann Chim</u>. 2003, Jan-Feb; 93(1-2): 1-10.
- Amstutz, H.C., Campbell, P., McKellop, H., Schmalzreid, T.P., Gillespie, W.J., Howie, D., Jacobs, J., Medley, J., Merritt, K.. "Metal on Metal Total Hip Replacement Workshop Consensus Document". Clin Orthop. 1996, Aug(329 Suppl): S297-303
- 3. Amstutz, H.C., Grigoris, P.. "*Metal on Metal Bearings in Hip Arthroplasty*". Clin Orthop. (329 Suppl): S11-34, Aug 1996.
- 4. Andrew, T.A., Berridge, D., Thomas, A., Duke, R.N.. "Long-term Review of Ring Total Hip Arthroplasty". Clin Orthop. (201): 111-122, Dec 1985.
- 5. Anissian, H.L., Stark, A., Gustafson, A., Good, V., Clarke, I.C.. "Metal-on-Metal Bearing in Hip Prosthesis Generates 100-Fold Less Wear Debris than Metal-on-Polyethylene". Acta Orthop Scand. 1999 Dec;70(6):578-82
- Archibeck, M.J., Jacobs, J.J., Black, J.. "Alternate Bearing Surfaces in Total Joint Arthorplasty: Biologic Considerations". Clin Orthop. 379:12-21, 2000.
- 7. August, A.C., Aldam, C.H., Pynsent, P.B.. "The McKee-Farrar Hip Arthroplasty. A Long-term Study". JBJS (BR). 86 (4): 520-527, Aug 1986.
- 8. Beaule, P.E., Dorey, F.J., LeDuff, M., Gruen, T., Amstutz, H.C.. "Risk Factors Affecting Outcome of Metal-on-Metal Surface Arthroplasty of the Hip". Clin Orthop. 2004 Jan; (418):87-93.
- 9. Benson, M.K., Goodwin, P.G., Brostoff, J.. "Metal Sensitivity in Patients with Joint Replacement Arthroplasties". BR Med J 4. (5993): 374-375, Nov 1975.
- 10. Besong, A.A., Lee, R., Farrar, R., Jin, Z.M.. "Contact Mechanics of a Novel Metal-on-Metal Total Hip Replacement", Proc Inst Mech Eng [H]. 2001;215(6):543-8.
- 11. Black, J.. "Metal on Metal Bearings. A Practical Alternative to Metal on Polyethylene Total Joints?". Clin Orthop. 1996 Aug;(329 Suppl):S244-55.
- 12. Boehler, N.. "Experiences with Metal on Metal Components in THR". Acta Orthop Belg. 1997;63 Suppl 1:96-7
- 13. Bohler, M., Kanz, F., Schwarz, B., Steffan, I., Walter, A., Plenk, H. Jr., Knahr, K.. "Adverse Tissue Reactions to Wear Particles from Co-Alloy Articulations, Increased by Alumina-Blasting Particle Contamination from Cementless Ti-based Total Hip Implants. A Report of Seven Revisions with Early Failure". JBJS Br. 2002 Jan;84(1):128-36.
- 14. Bohler, N.. "Metal/Metal Articulating Interfaces". Orthopedics. 1995 Sep;18(9):879-80.

- 15. Brodner, W., Bitzan, P., Meisinger, V., Kaider, A., Gottsauner-Wolf, F, Kotz, R.. "Serum Cobalt Levels After Metal-on-Metal Total Hip Arthroplasty". JBJS Am. 2003 Nov;85-A(11):2168-73.
- 16. Brodner, W., Bitzan, P., Meisinger, V., Kaider, A., Gottsauner-Wolf, F., Kotz, R.. "Elevated Serum Cobalt with Metal-on-Metal Articulating Surfaces". JBJS BR. 79 (5): 885, Sep 1997.
- 17. Brodner, W., Grohs, J.G., Bitzan, P., Meisinger, V., Kovarik, J., Kotz, R.. "Serum Cobalt and Serum Chromium Level in 2 Patients with Chronic Renal Failure after Total Hip Prosthesis Implantation with Metal-Metal Gliding Contact". Z Orthop Ihre Grenzgeb. 2000 Sep-Oct;138(5):425-9.
- 18. Browne, M., Gregson, P.J.. "Effect of Mechanical Surface Pretreatment in Metal Ion Release". Biomaterials. 2000 Feb;21(4):385-92.
- 19. Bryant, M.J., Mollan, R.A., Nixon, J.R.. "Survivorship Analysis of the Ring Hip Arthroplasty". J. Arthoplasty. (6 Suppl): S5-10, 1991.
- 20. Campbell, P.A., Wang, M., Amstutz, H.C., Goodman, S.B.. "Positive Cytokine Production in Failed Metal-on-Metal Total Hip Replacements". Acta Orthop Scand. 2002 Oct;73(5):506-12.
- 21. Chan, F.W., Bobyn, J.D., Medley, J.B., Krygier, J.J., Yue, S., Tanzer, M.. "Engineering Issues and Wear Performances of Metal on Metal Hip Implants". Clin Orthop. (333): 96-107, Dec 1996.
- 22. Clarke, M.T., Lee, P.T., Arora, A., Villar, R.N.. "Levels of Metal Ions After Small- and Large-Diameter Metal-on-Metal Hip Arthroplasty". JBJS Br. 2003 Aug;85(6):913-7.
- 23. Clarke, M.T., Lee, P.T., Villar, R.N.. "Dislocation After Total Hip Replacement in Relation to Metal-on-Metal Bearing Surfaces". JBJS Br. 2003 Jul;85(5);650-4.
- 24. Coleman, M.P.. "Cancer Risk from Orthopedic Prostheses". Ann Clin Lab Sci. 1996 Mar-Apr;26(2):139-46.
- 25. Daniel, J., Pynsent, P.B., McMinn D.J.. "Metal-on-Metal Resurfacing of the Hip in Patients Under the Age of 55 Years with Osteoarthritis". JBJS Br. 2004 Mar;86(2)177-84.
- 26. Djerf, K., Wahlstrom, O.. "Total Hip Replacement Comparison Between the McKee-Farrar and Charnley Prostheses in a 5-year Follow-up Study". Arch Orthop Trauma Surg. 105 (3): 158-162, 1986.
- 27. Doorn, P.F., Campbell, P.A., Amstutz, H.C.. "Metal versus Polyethylene Wear Particles in Total Hip Replacements. A Review". Clin Orthop. (329 Suppl): S206-216, Aug 1996.
- 28. Doorn, P.F., Campbell, P.A., Worrall, J., Benya, P.D., McKellop, H.A., Amstutz, H.C.. "Metal Wear Particle Characterization from Metal on Metal Total Hip Replacements: Transmission Electron Microscopy Study of Periprosthetic Tissues and Isolated Particles". J Biomet Mater Res. 42 (1): 103-111, Oct 1998.
- 29. Doorn, P.F., Mirra J.M., Campbell, P.A., Amstutz, H.C.. "Tissue Reaction to Metal on Metal Hip Prostheses". Clin Orthop. 1996 Aug;(329 Suppl):S187-205.

- 30. Dorr, L.D., Hilton, K.R., Wan, Z., Markovich, G.D., Bloebaum, R.. "Modern Metal on Metal Articulation for Total Hip Replacements". Clin Orthop. (333): 108-117, Dec 1996.
- 31. Dorr, L.D., Wan, Z., Longjohn, D.B., Dibois, B., Murken, R.. "Total Hip Arthroplasty with Use of the Metasul Metal-on-Metal Articulation. Four to Seven-Year Results". JBJS Am. 2000 Jun;82(6):789-98.
- 32. Evans, E.M.. "Metal Sensitivity as a Cause of Bone Necrosis and Loosening of the Prosthesis in Total Joint Replacement". JBJS 56-B.(4):626-642, 1974.
- 33. Gillespie, W.J., Frampton, C.M., Henderson, R.J., Ryan, P.M.. "The Incidence of Cancer Following Total Hip Replacement". <u>JBJS 70-B</u>. (4):539-542, 1988.
- 34. Gleizes, V., Poupon, J., Lazennec, J.Y., Chamberlin, B., Saillant, G.. "Value and Limits of Determining Serum Cobalt Levels in Patients with Metal on Metal Articulating Prostheses". Rev Chir Orthop Reparatrice Appar Mot. 85 (3): 217-225, Jun 1999.
- 35. Goldsmith, A.A., Dowson, D., Isaac, G.H., Lancaster, J.G.. "A Comparative Joint Simulator Study of the Wear of Metal-on-Metal and Alternative Material Combination in Hip Replacements". Proc Inst Mech Eng [H]. 2000;214(1):39-47.
- 36. Hallab, N., Merritt, K., Jacobs, J.J.. "Metal Sensitivity in Patients with Orthopaedic Implants". JBJS 83-A. (3):428-436, 2001.
- 37. Hallab, N.J., Messina, C., Skipor, A., Jacobs J.J.. "Differences in the Fretting Corrosion of Metal-Metal Modular Junctions of Total Hip Replacements". J Orthop Res. 2004 Mar;22(2):250-9.
- 38. Hallab, N.J., Mikecz, K., Jacobs, J.J.. "A Triple Assay Technique for the Evaluation of Metal-induced Delayed-type Hypersensitivity Responses in Patients With or Receiving Total Joint Arthroplasty". J Biomed Mater Res. 53(5):480-489, 2000.
- 39. Hallab, N.J., Mikecz, K., Vermes, C., Skipor, A., Jacobs, J.J.. "Orthopaedic Implant Related Metal Toxicity in Terms of Human Lymphocyte Reactivity to Metal-Protein Complexes Produced from Cobalt-base and Titanium-base Implant Alloy Degradation". Mol Cell Biochem. 222(1-2):127-136, 2001.
- 40. Hallab, N.J., Mikecz, K., Vermes, C., Skipor, A., Jacobs, J.J.. "Differential Lymphocytic Reactivity to Serum-derived Metal-Protein Complexes Produced from Cobalt-based and Titanium-based Implant Alloy Degradation" J Biomed Mater Res. 56(3):427-436, 2001.
- 41. Hallab, N.J., Vermes, C., Messina, C., Roebuck, K.A., Glant, T.T., Jacobs, J.J.. "Concentration- and Composition-dependent Effects of Metal Ions on Human MG-63 Osteoblasts". J Biomed Mater Res. 60(3):420-433, 2002.
- 42. Hallam, P., Haddad, F., Cobb, J.. "Pain in the Well-Fixed, Aseptic Titanium Hip Replacement. The Role of Corrosion". JBJS Br. 2004 Jan;86(1):27-30.
- 43. Harding, I., Bonomo, A., Crawford, R., Psychoyios, V., Delves, T., Murray, D., McLardy-Smith, P.. "Serum Levels of Cobalt and Chromium in a Complex Modular Total Hip Arthroplasty System". J Arthroplasty. 2002 Oct;17(7):893-5.

- 44. Hennig, F.F., Raithel, H.J., Schaller, K.H., Dohler, J.R.. "Nickel-, Chrom- and Cobalt-Concentrations in Human Tissue and Body Fluids of Hip Prosthesis Patients". J Trace Elem Electrolytes Health Dis. 1992 Dec;6(4):239-43.
- 45. Higuchi, F., Inoue, A., Semlitsch, M.. "Metal-on-Metal CoCrMo McKee-Farrar Total Hip Arthroplasty: Characteristics from a Long-term Follow-up Study". Arch Orthop Trauma Surg. 116(3): 121-4, 1997.
- 46. Hilton, K.R., Dorr, L.D., Wan, Z., McPherson, E.J.. "Contemporary Total Hip Replacement with Metal on Metal Articulation". Clin Orthop. 1996 Aug;(329Suppl):S99-105.
- 47. Howie, D.W., Rogers, S.D., McGee, M.A., Haynes, D.R.. "Biologic Effects of Cobalt Chrome in Cell and Animal Models". Clin Orthop. 329 Suppl:S217-232, 1996.
- 48. Inzerillo, V.C., Garino, J.P.. "Alternative Bearing Surfaces in Total Hip Arthroplasty". J South Orthop Assoc. 2003 Summer;12(2):106-11.
- 49. Jacobs, J.J., Skipor, A.K., Black, J. Urban, R., Galante, J.O.. "Release and Excretion of Metal in Patients Who Have a Total Hip Replacement Component Made of Titanium-base Alloy". JBJS Am. 1991 Dec;73(10):1475-86.
- 50. Jacobs, J.J., Skipor, A.K., Doorn, P.F., Campbell, P., Schmalzried, T.P., Black, J. and Amstutz, H.C.. "Cobalt and Chromium Concentrations in Patients with Metal on Metal Total Hip Replacements". Clin Orthop. (329 Suppl): S256-263, Aug 1996.
- 51. Jacobs, J.J., Skipor, A.K., Patterson, L.M., Hallab, N.J., Paprosky, W.G., Black, J., Galante, J.O.. "Metal Release in Patients Who Have Had a Primary Total Hip Arthroplasty. A Prospective, Controlled, Longitudinal Study". JBJS Am. 1998 Oct;80(10):1447-58.
- 52. Jacobsson, S.A., Djerf, K., Wahlstrom, O.. "Twenty-year Results of McKee-Farrar versus Charnley Prosthesis". Clin Orthop. (329 Suppl): S60-68, Aug 1996.
- 53. Jazrawi, L.M., Kummer, F.J., Di Cesare, P.E.. "Hard Bearing Surfaces in Total Hip Arthroplasty". Am J Orthop. 27 (4): 283-92, Apr 1998.
- 54. Jazrawi, L.M., Kummer, F.J., DiCesare, P.E.. "Alternative Bearing Surfaces for Total Joint Arthroplasty". J Am Acad Orthop Surg. 1998 Jul-Aug;6(4):198-203.
- 55. Korovessis, P., Petsinis, G., Repanti, M.. "Zweymueller with Metal-on-Metal Articulation: Clinical Radiological and Histological Analysis of Short-Term Results". Arch Orthop Trauma Surg. 2003 Feb; 123(1):5-11. Epub 2002 Dec 19.
- 56. Kossovsky, N., Liao, K., Millett, D., Feng, D., Campbell, P.A., Amstutz, H.C., Finerman, G.A., Thomas, B.J., Kilgus, D.J, Cracchiolo, A., et all.. "Periprosthetic Chronic Inflammation Characterized through the Measurement of Superoxide Anion Production by Synovial-derived Macrophages" Clin Orthop. 263:263-71, 1991.
- 57. Kreibich, D.N., Moran, C.G., Delves, H.T., Owen, T.D., Pinder, I.M.. "Systemic Release of Cobalt and Chromium After Uncernented Total Hip Replacement". JBJS Br. 1996
  Jan;78(1):18-21.

- 58. Langkamer, V.G., Case, C.P., Heap, P., Taylor, A., Collins, C., Pearse, M., Solomon, L.. "Systemic Distribution of Wear Debris After Hip Replacement. A Cause for Concern?". JBJS 74-B. (6):831-839, 1992.
- 59. Lee, S.H., Brennan, F.R., Jacobs, J.J., Urban, R.M., Ragasa, D.R., Glant, T.T.. "Human Monocyte/Macrophages Response to Cobalt-Chromium Corrosion Products and Titanium Particles in Patients with Total Joint Replacement". J Orthop Res. 1997 Jan; 15(1):40-9.
- 60. Lhotka, C., Szekeres, T., Steffan, I., Zhuber, K., Zweymuller, K.. "Four-Year Study of Cobalt and Chromium Blood Levels in Patients Managed with Two Different Metal-on-Metal Total Hip Replacements". <u>J Orthop Res</u>. 2003 Mar;21(2):189-95.
- 61. Lindholm, T.S., Viljakka, T., Vankka, E., Popov, L., Lindholm, T.C.. "Development of Heterotopic Ossification Around the Hip. A Long-term Follow-up of Patients Who Underwent Surgery with Two Different Types of Endoprostheses". Arch Orthop Trauma Surg. 105 (5): 236-267, 1986.
- Lombardi, A.V. Jr., Mallory, T.H., Alexiades, M.M., Cuckler, J.M., Faris, P.M., Jaffe, K.A., Keating, E.M., Nelson, C.L. Jr., Ranawat, C.S., Williams, J., Wixson, R., Hartman, J.F., Capps, S.G., Kefauver, C.A.. "Short-Term Results of the M2a-Taper Metal-on-Metal Articulation". J Arthroplasty. 2001 Dec;16(8 Suppl 1): 122-8.
- 63. Lombardi, A.V. Jr., Mallory, T.H., Alexiades, M.M., Cuckler, J.M., Faris, P.M., Jaffe, K.A., Keating, E.M., Nelson, C.L. Jr., Ranawat, C.S., Williams, J., Wixson, R., Berend, K.R., Dodds, K.L., Adams, J.B.. "Mid-term Results of a Polyethylene-free Metal on Metal Articulation". J Arthroplasty. 2004 Feb;1992:260
- 64. MacDonald, S.J., McCalden, R.W., Chess, D.G., Bourne, R.B., Rorabeck, C.H., Cleland, D., Leung, F.. "Metal-on-Metal Versus Polyethylene in Hip Arthroplasty: A Randomized Clinical Trial". Clin Orthop. 2003 Jan; (406):282-96.
- 65. Mathiesen, E.B., Ahlbom, A., Bermann, F., Lingren, J.U.. "Total Hip Replacement and Cancer. A Cohort Study". JBJS 77-B. (3):345-350, 1995.
- McKellop, H., Park, S.H., Chiesa, R., Doorn, P., Lu, B., Normand, P., Grigoris, P., Amstutz, H.. "IN VIVO Wear of Three Types of Metal on Metal Hip Prostheses During Two Decades of Use". Clin Orthop. (329 Suppl): S128-140, Aug 1996.
- 67. McMinn, D., Treacy, R., Lin, K., Pynsent, P.. "Metal on Metal Surface Replacement of the Hip. Experience of the McMinn Prosthesis". Clin Orthop. (329 Suppl): S128-140, Aug 1996.
- 68. Merritt, K., Brown, S.A.. "Distribution of Cobalt-Chromium Wear and Corrosion Products and Biological Reactions". Clin Orthop. 1996 Aug ;(329 Suppl):S233-43.
- 69. Merritt, K., Rodrigo, J.J.. "Immune Response to Synthetic Materials. Sensitization of Patients Receiving Orthopaedics Implants". Clin Orthop. 326:71-79, 1996.
- 70. Nyren, O., McLaughlin, J.K., Gridley, G., Ekbom, A., Johnell, O., Fraumeni, J.F. Jr., Adami, H.O.. "Cancer Risk After Hip Replacement with Metal Implants: A Population-Based Cohort Study in Sweden". J Natl Cancer Inst. 1995 Jan 4;87(1):28-33.


- 71. Olsen, R.V., Munk, P.L., Lee, M.J., Janzen, D.L. MacKay, A.L., Xiang, Q.S., Marsi, B. "Metal Artifact Reduction Sequence: Early Clinical Applications". Radiographics. 2000 May-Jun;20(3):699-712.
- 72. Paavolainen, P., Pukkala, E., Pulkkinen, P., Visuri, T.. "Cancer Incidence in Finnish Hip Replacement Patients from 1980 to 1995: A Nationwide Cohort Study Involving 31,651 Patients". J Arthroplasty. 14(3):272-280, 1999.
- 73. Pazzaglia, U.E., Minoia, C., Gualtieri, G., Gualtieri, I., Riccardi, C., Ceciliani, L.. "Metal Ions in Body Fluids After Arthroplasty". Acta Orthop Scand. 1986 Oct;57(5):415-8.
- 74. Randle, R. and Goriev, K.. "Metal-on-Metal Articulation in Total Hip Arthroplasty Preliminary Results in 57 Cases". Aust N Z J Surg. 67 (9): 634-636, Sept 1997.
- 75. Ray, A.. "15 Years Survival of the Mac-Kee Farrar Metal Hip Prosthesis. Apropos of 58 Cases and 4 Explanted Cups". Rev Chir Orthop Reparatrice Appar Mot. 82 (1): 85-89, 1996.
- 76. Reiker, C., Konrad, R., Schon, R.. "IN VITRO Comparison of the Two Hard-Hard Articulations for Total Hip Replacement". Proc Inst Mech Eng [H]. 2001;215(2):153-60.
- 77. Reinisch, G., Judmann, K.P., Lhotka, C., Lintner, F., Zweymuller, K.A.. "Retrieval Study of Uncemented Metal-Metal Hip Prostheses Revised for Early Loosening". Biomaterials. 2003 Mar;24(6):1081-91.
- 78. Rotem, A.. "Effect of Implant Material Properties on the Performance of a Hip Joint Replacement". J Med Eng Technol. 1994 Nov-Dec; 18(6):208-17.
- 79. Santavirta, S., Bohler, M., Harris, W.H., Konttinen, Y.T., Lappalainen, R., Muratoglu, O., Rieker, C., Salzer, M.. "Alternative Materials to Improve Total Hip Replacements Tribiology". Acta Orthop Scand. 2003 Aug;75(4):380-8.
- 80. Savarino, L., Granchi, D., Ciapetti, G., Cenni, E., Greco, M., Rotini, R., Veronesi, C.A., Baldini, N., Giunti, A.. "Ion Release in Stable Hip Arthroplasties Using Metal-on-Metal Articulating Surfaces: A Comparison Between Short- and Medium-Term Results".

  J Biomed Mater Res. 2003 Sep 1;66A(3):450-6.
- 81. Savarino, L., Granchi, D., Ciapetti, G., Cenni, E., Nardi Pantoli, A., Rotini, R., Veronesi, C.A., Baldini, N., Giunti, A.. "Ion Release in Patients with Metal-on-Metal Hip Bearings in Total Joint Replacement: A Comparison with Metal-on-Polyethylene Bearings". J Biomed Mater Res. 2002;63(5):467-74.
- 82. Schaffer, A.W., Pilger, A., Engelhardt, C., Zweymueller, K., Ruediger, H.W.. "Increased Blood Cobalt and Chromium After Total Hip Replacement". <u>J Toxicol Clin Toxicol</u>. 1999;37(7):839-44.
- 83. Schmalzried, T.P, Szuszczewicz, E.S., Akizuki, K.H., Petersen, T.D., Amstutz, H.C.. "Factors Correlating with Long Term Survival of McKee-Farrar Total Hip Prostheses". Clin Orthop. (329 Suppl): S48-59, Aug 1996.
- 84. Schmalzried, T.P., Fowble, V.A, Ure, K.J., Amstutz, H.C.. "Metal on Metal Surface Replacement of the Hip. Technique, Fixation, and Early Results". Clin Orthop. 1996 Aug;(329 Suppl):S106-14.

- 85. Schmalzried, T.P., Peters, P.C., Maurer, B.T., Bragdon, C.R., Harris, W.H.. "Long-duration Metal-on-Metal Total Hip Arthroplasties with Low Wear of the Articulation Surfaces". J. Arthroplasty. 11 (3): 322-331, Apr 1996.
- 86. Sieber, H.P., Rieker, C.B., Kottig, P.. "Analysis of 118 Second Generation Metal-on-Metal Retrieved Hip Implants". JBJS Br. 1999 Jan;81(1):46-50.
- 87. Smith, S.L., Dowson, D., Goldsmith, A.A.. "The Effect of Femoral Head Diameter Upon Lubrication and Wear of Metal-on-Metal Total Hip Replacements". Proc Inst Mech Eng [H]. 2001;215(2):161-70.
- 88. St. John, K.R., Zardiackas, L.D., Poggie, R.A.. "Wear Evaluation of Cobalt-Chromium Alloy for Use in a Metal-on-Metal Hip Prosthesis". J Biomed Mater Res. 2004 Jan 15;68B(1):1-14.
- 89. Streicher, R.M.. "Metal-on-Metal Articulation in Total Hip Arthroplasty: The Case for Using Metal-on-Metal". <u>J Arthroplasty</u>. 13(3): 343-5, Apr 1998.
- 90. Szuszczewicz, E.S., Schmalzried, T.P. and Petersen, T.D.. "*Progressive Bilateral Pelvic Osteolysis in a Patient with McKee-Farrar Metal-Metal Total Hip Prosthesis*".

  <u>J Arthroplasty</u>. 12 (7): 819-824, Oct 1997.
- 91. Tager, G., Euler, E., Plitz, W.. "Changes in Shape of the McKee-Farrar Hip Endoprosthesis". Orthopade. 26 (2): 142-151, Feb 1997.
- 92. Tharani, R., Dorey, J.F., Schmalzried, T.P.. "The Risk of Cancer Following Total Hip or Knee Arthroplasty". JBJS 83-A. (5):774-780, 2001.
- 93. Vale, L., Wyness, L., McCormack, K., McKenzie, L., Brazzelli, M., Stearns, S.C.. "A Systemic Review of the Effectiveness and Cost-Effectiveness of Metal-on-Metal Hip Resurfacing Arthroplasty for Treatment of Hip Disease". Health Technol Assess. 2002;6(15):1-109.
- 94. Visuri, T. and Koskenvuo, M.. "Cancer Risk After McKee-Farrar Total Hip Replacement". Orthopedics. 14 (2): 137-142, Feb 1991.
- 95. Visuri, T., Pukkala, E., Paavolainen, P., Pulkkinen, P. and Riska, E.B.. "Cancer Risk After Metal on Metal and Polyethylene on Metal Total Hip Arthroplasty". Clin Orthop. (329 Suppl): S280-289, Aug 1996
- 96. Wagner, M., Wagner, H.. "Medium-Term Results of a Modern Metal-on-Metal System in Total Hip Replacement". Clin Orthop. 2000 Oct;(379):123-33.
- 97. Wagner, M., Wagner, H.. "Preliminary Results of Uncemented Metal on Metal Stemmed and Resurfacing Hip Replacement Arthroplasty". Clin Orthop. 1996 Aug;(329 Suppl):S78-88.
- 98. Weber, B.G., Fiechter, T.. "Polyethylene Wear and Late Loosening of a Total Prosthesis of the Hip Joint. New Perspectives for Metal/Metal Pairing of the Capsule". Orthopade. 18 (5): 370-376, Sept 1989.

- 99. Weber, B.G., Semlitsch, M.F., Streicher, R.M.. "Total Hip Joint Replacement Using a CoCrMo Metal-Metal Sliding Pairing". Nippon Seikeigeka Gakkai Zasshi. 1993 May;67(5):391-8.
- 100. Weber, B.G.. "Experience with the Metasul Total Hip Bearing System". Clin Orthop. (329 Suppl): S69-77, Aug 1996.
- Willert, H.G., Buchhorn, G.H., Gobel, D., Koster, G., Schaffner, S., Schenk, R., Semlitsch, M.. "Wear Behavior and Histopathology of Classic Cemented Metal on Metal Hip Endoprostheses". Clin Orthop. (329 Suppl): S160-186, Aug 1996.
- 102. Willert, H.G., Semlitsch, M.. "Tissue Reactions to Plastic and Metallic Wear Products of Joint Endoprostheses". Clin Orthop. 333:4-14, 1996.
- 103. Zaoussis, A.L. and Patikas, A.F.. "Experience with Total Hip Arthroplasty in Greece, The First 20 Years. A Particular Reference to Long-term Results with the McKee-Farrar Technique". Clin Orthop. 246: 39-47, Sept 1989.



Entrez PubMed

Overview
Help | FAQ
Tutorial
New/Noteworthy
E-Utilities

PubMed Services
Journals Database
MeSH Database
Single Citation Matcher
Batch Citation Matcher
Clinical Queries
LinkOut
Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

Cobalt blood levels after total hip replacement (THR): a new follow-up study in Trieste (Italy).

Adami G, Smarrelli D, Martinelli B, Acquavita A, Reisenhofer E.

Department of Chemical Sciences, University of Trieste, Via Giorgieri 1, 34127 Trieste, Italy. adami@dsch.univ.trieste.it

Cobalt and chromium release in patients undergone a metal-on-metal total hip replacement (THR) is a matter recently discussed and whose we do not have enough information about it yet. In literature there is little data and not completely indicative, in the biological fluids and organs the amounts of released metals are different. This is also due to the fact that cobalt and chromium blood levels can change depending on physical and working activity, individual feeding and metabolism. The results obtained confirm the presence of an increase of cobalt inthe blood of patients after total hip replacement, while the chromium levels are almost alike: average values in patients operated are 4.1 +/- 1.5 microg/L for cobalt (0.3 +/- 0.1 microg/L in the control group) and 4.5 +/- 2.9 microg/L for chromium (4.7 +/- 2.4 microg/L in the control group). In spite of the cobalt values stand below the concentration generally considered dangerous, the difference between the two examined groups points out that a risk exists for the health of these patients. These results must be confirmed by further studies, providing better information and more reliable and biocompatible materials.

PMID: 12650568 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort Send to Text

Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer

Apr 19 2004 06:53:45


**NLM Gateway** 

Clinical Alerts
ClinicalTrials.gov

PubMed Central

Privacy Policy

TOXNET Consumer Health



PMID: 8769344 [PubMed - indexed for MEDLINE]

consensus statements were developed.

Display Abstract Show: 20 Sort Send to Text

clinical, tribologic, and biologic considerations of all metal bearings. Based on the

of the panelists, and the current regulatory, legal, and economic environment,

scientific presentations at the symposium, the extant literature, the clinical experience

Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer

Apr 19 2004 06:53:45







				•			****	Journals	Books
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	MiMO	PMC Go C	Journals Dear	5001111
Search	PubMed	for Limits	Previe	w/Index	History	С	lipboard	Deta	ails
About Enti	rez	<b>Display</b> Abstra	ect	Show	w: 20 ¥ Sor	t s	Send to	Text	2

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

**PubMed Services** Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents **NLM Gateway** TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

Related Articles, Links

Comment in:

Clin Orthop. 1997 Jul; (340):283-4.

## Metal on metal bearings in hip arthroplasty.

Amstutz HC, Grigoris P.

Joint Replacement Institute, Orthopaedic Hospital, Los Angeles, CA, USA.

Periprosthetic osteolysis caused by wear debris released from the bearing surface of polyethylene components is the major problem in contemporary hip arthroplasty. Several types of metal on metal prostheses were developed in the 1960s, but by the mid 1970s they were completely displaced by polyethylene bearings. There have been several generations of all metal components with significant variation in design, tolerances, and bearing surface quality. A number of these hips have survived for more than 25 years because of low wear rates and minimal osteolysis. Identification of the characteristics that contributed to long term function is important. The historical development and clinical results of metal on metal hip arthroplasties are presented. Factors that led to the abandonment of the metal on metal bearings are related to: (1) the early success of the Charnley prosthesis; (2) the frictional torque issue; (3) carcinogenesis concerns; (4) metal sensitivity concerns; (5) high infection rates; and (6) increased strain rates in periprosthetic bone and fatigue fractures of the acetabular floor. The accumulated experience to date enables one to evaluate all the factors with a different perspective and makes the use of newer metal on metal bearings a viable option in younger patients.

#### Publication Types:

- Review
- Review, Multicase

PMID: 8769320 [PubMed - indexed for MEDLINE]


<u></u>	
Display Abstract	Show: 20 Sort Send to Text

Write to the Help Desk NCBI NLM NIH Department of Health & Human Services Freedom of Information Act | Disclaimer

Jun 1 2004 06:47:16







			•	•		
Entrez	PubMed	Nucleotide	Protein Genon	ne Structure	OMIM PM	lood alamuol Book
Search	PubMed	for		· · · · · · · · · · · · · · · · · · ·	Go	Clear
Search	America de Companyo de Compaño de	Limits	Preview/Index	History	Clipboa	rd Details
About Entr	ez	<b>Display</b> Abstra	ct 🕌	Show: 20	Sort 👺 Se	end to Text

About Entrez

**Text Version** 

Entrez PubMed Overview Help I FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents **NLM Gateway** TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ 1: Clin Orthop. 1985 Dec;(201):111-22.

Related Articles, Links

## Long-term review of ring total hip arthroplasty.

Andrew TA, Berridge D, Thomas A, Duke RN.


In a five- to 12-year follow-up survey of 179 sequential cementless Ring metal-on-metal total hip arthroplasties, 55 (31%) of the patients were found to have died as a result of nonorthopedic conditions. Analysis of the records demonstrated that 20% of these patients had had poor results attributable to pain. Of the remaining 124 patients, 116 (94%) attended for full clinical and radiologic review yielding a total of 154 hips. Using Ring's classification, 75 hips were judged to have excellent or good results. Forty-one hips were graded as fair or poor as a result of pain, and an additional 15 hips were revised for symptomatic loosening. There were five cases of Brooker Grade IV periarticular ossification, four cases of gross metal reaction requiring prosthetic removal, and two cases of infection. There was considerable variation in the radiographic appearance of the hips, and at times radiographic changes were inconsistent with clinical symptoms. Eleven of the revised hips were converted to longer and largerdiameter uncemented Ring femoral components. Nine of these yielded only fair or poor results at the time of review, whereas both cases in which the femoral component was cemented were associated with good results.

PMID: 4064395 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort Send to Text

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer

Jun 1 2004 06:47:16



Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services
Journals Database
MeSH Database
Single Citation Matcher
Batch Citation Matcher
Clinical Queries
LinkOut
Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

Metal-on-metal bearing in hip prosthesis generates 100-fold less wear debris than metal-on-polyethylene.

Anissian HL, Stark A, Gustafson A, Good V, Clarke IC.

Department of Orthopaedic Surgery, Karolinska Hospital, Stockholm, Sweden. Anissian@ort.ks.se

Aseptic loosening due to osteolysis in total hip replacement has been related to wear debris released from prosthetic components. Retrospective longterm observations of patients with the metal-on-metal prosthesis has shown long-term survivorship and good mechanical performance. Thus, the new and modified metal-on-metal prosthesis has been introduced on the market. Historical clinical data from the 1st generation metal-on-metal hip prosthesis may not be relevant for the 2nd generation of metal-on-metal hip prosthesis. Therefore, preclinical testing of the prosthesis must be conducted before clinical evaluation. We assessed the tribological performance of the metal-on-metal prosthesis versus the metal-on-polyethylene prosthesis introduced on the market as Metasul and Protasul, respectively. In a 12-channel joint simulator, 6 metal-on-metal bearing and 3 metal on polyethylene prostheses were tested, with the same number of corresponding soak controls. The wear was assessed gravimetrically. The "steady-state" wear-rates from the metal-on-metal prosthesis were almost 100 times less than that from the metal-on-polyethylene prosthesis. The tribological wear performance of the metal-on-metal hip prosthetic system is promising.


PMID: 10665722 [PubMed - indexed for MEDLINE]

<u> </u>		
Display Abstract	Show: 20 Sort Sort	end to Text


Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer

Apr 19 2004 06:53:45

OCCETTO FINITION







Send to

Entrez	PubMed	Nucleotide	Protein	Genome	Structure	MIMO	PMC	Journals	В
Search	PubMed	<b>∃</b> for	n a contragación (o assessas C desimbles processoradas)	alle for the first of the second construction of the foreign construction of the second construction o	and seeds, because the measure of the time and as 3 th, responsible daily through	and the second section of the section of	Go] Cl	ear	talla
	E-recomp companies () () with interest	Limits	Previ	ew/index	History	Cli	pboard	De	tails
							white test, who should	4 th	

Show: 20

About Entrez

**Text Version** 

E-Utilities

Cubby

Abstract

Display

Related Articles, Links

Text

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy

PubMed Services
Journals Database
MeSH Database
Single Citation Matcher
Batch Citation Matcher
Clinical Queries
LinkOut

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ 1: Clin Orthop. 2000 Oct;(379):12-21.

Alternate bearing surfaces in total joint arthroplasty: biologic considerations.

Sort

Archibeck MJ, Jacobs JJ, Black J.

Department of Orthopaedic Surgery, Rush-Presbyterian-St Luke's Medical Center, Chicago, IL 60612, USA.

The problem of periprosthetic osteolysis is currently the major limiting factor in joint arthroplasty longevity. Because this process has been shown to be primarily a biologic response to wear particles, corrosion products, or both, efforts to reduce particle generation are being undertaken. These efforts include the development of modified polyethylene and alternative articulating surfaces. These alternate bearing surfaces currently include ceramic-on-polyethylene, ceramic-on-ceramic, and metal-on-metal. Although these alternate bearings diminish or eliminate the generation of polyethylene particles, ceramic and metal particles are produced. The purpose of the current review is to discuss the literature that addresses the biologic response to these particles, locally and systemically.

Publication Types:

- Review
- Review, Academic

PMID: 11039787 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort	Send to Te	ext 🚬 🛨

Write to the Help Desk

NCBI | NLM | NIH

Department of Health & Human Services

Freedom of Information Act | Disclaimer

Jun 1 2004 06:47:16

Privacy Policy







8	NCBI		Put	OMe	d		ibrary dicine		
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	OMIMO	PMC	Journals	Books
Search	PubMed	for Limits	Pre	view/Index	History	С	Go Cl	<b>ear</b> Det	ails
About Entre	ZZ	Display Ab	stract	Sho	w: 20 🔁 Sc	ort 🔣	Send to	Text	
Entrez Pu Overview Help   FAQ Tutorial		□ <b>1:</b> J Bone J			8(4):520-7. <b>throplasty</b>	. A long	-term s	•	rticles, Links
New/Notew E-Utilities	rorthy	August	AC, Aldai	n CH, Pynser	nt PB.				
	atabase abase tion Matcher tion Matcher	arthropl survivin were go .score w was und	asties were g arthropla ood or exce ith 78% of dertaken ar	e performed in sties have bee llent results in	of 808 McKee-f the Norfolk ar in reviewed at 49% of the ar ittle or no pair ip study of 81	nd Norwick average f throplasti n. A compl	n Hospitai. ollow-up o es as judg rehensive	of these, 2 of 13.9 years ed by the H radiographic	. There arris hip analysis
Related F Order Doc NLM Gate TOXNET Consumer	way	present PMID: 3		ubMed - index	ed for MEDLIN	IE]			
Clinical Ale ClinicalTria PubMed C	erts als.gov	Display At	ostract	Sho	<sub>ow</sub> . 20 😿 S	ort 🗽	Send t	o Text	

Write to the Help Desk


NCBI | NLM | NIH

Department of Health & Human Services

Freedom of Information Act | Disclaimer







Mary San San San			7	•		
Entrez	PubMed	Nucleotide	Protein Genor	rie Structure	OMIM PMC	Journals Books
Search F	ubMed	for			I Go	Clear
Jeanen ,		Limits	Preview/Index	History	Clipboard	Details
About Entrez		Display Abstra	ct 🚆	Show: 20 S	ort 🔻 Send	to Text

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services
Journals Database
MeSH Database
Single Citation Matcher
Batch Citation Matcher
Clinical Queries
LinkOut
Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ **1**: Clin Orthop. 2004 Jan;(418):87-93.

Related Articles, Links



## Risk factors affecting outcome of metal-on-metal surface arthroplasty of the hip.

Beaule PE, Dorey FJ, LeDuff M, Gruen T, Amstutz HC.

David Geffen School of Medicine at UCLA, Joint Replacement Institute at Orthopaedic Hospital, Los Angeles, CA 90007, USA. pbeaule@mednet.ucla.edu

We evaluated radiologic and clinical features affecting the outcome of hybrid metal-onmetal surface arthroplasty of the hip in 119 hips in patients 40 years and younger. Only the hips that had either failed or had minimum 2-year followup were reviewed. Ninetyfour hips in 83 patients with a mean age of 34.2 years (range, 15-40 years) were reviewed. Seventy-one percent of the patients were males and 29% of the patients were females; 14% had previous surgery. The Chandler index and surface arthroplasty risk index were calculated. The mean followup at 3 years (range, 2-5 years) showed that three hips were converted to a total hip replacement at a mean of 27 months (range, 2-50 months) after the original surgery, and 10 hips had significant radiologic changes. The mean surface arthroplasty risk index for these 13 problematic hips versus the remaining hips was significantly higher, 4.7 and 2.6, respectively. The mean angle between the prosthesis stem and femoral shaft in the problematic group was significantly smaller than in the remaining hips: 133 degrees and 139 degrees, respectively. With a surface arthroplasty risk index score greater than 3 the relative risk of early problems is 12 times greater than if surface arthroplasty risk index less than or equal to 3.

PMID: 15043098 [PubMed - in process]

	<del></del>
Display Abstract	Show: 20 ► Sort ► Send to Text ►

Write to the Help Desk
NCBI | NLM | NiH
Department of Health & Human Services
Freedom of Information Act | Disclaimer

Privacy Policy







8	NCBI		Pub	Me	ed	Libra of Medic		
Entrez	PubMed	Nuclectice	Protein	Genome	Structure	- Section	MC Journ	nais Books
Search F	PubMed	for Limits	Prev	view/Index	History	Clipbo	ming specialism and see	Details
About Entrez	:	<b>Display</b> Ab	stract	Sho	w: 20 💟 So	rt 🕎 🎉	Send to Te	xt 🔄
Text Version								
Entrez Pub Overview Help   FAQ Tutorial New/Notewo E-Utilities				ity in patio	74-5. ents with jo	oint replac	•	elated Articles, Links
PubMed So Journals Da MeSH Datal Single Citati Batch Citatio Clinical Que LinkOut Cubby	tabase base on Matcher on Matcher	A high i metal (l	ncidence of McKee) hip	arthroplasties dence of met	netal sensitivity The Patients with The Sensitivity the The Patient of the pre-	metal-to-plas an a control g	roup awaitir	ng operation. If
Related Re Order Docu NLM Gatew TOXNET Consumer I Clinical Alei ClinicalTrial PubMed Ce	ments ray Health rts s.gov		1192078 [Pu	ubMed - index	ed for MEDLIN		Send to Te	ext


Write to the Help Desk

NCBI | NLM | NIH

Department of Health & Human Services
Freedom of Information Act | Disclaimer







Entrez	PubMed	Nucleotide	Protein	Genome	Structure	OMIM	PMC 5 September 1 September 1	Journals Posteron	Book
Tr.	PubMed	for					Go C	lear	
Search L		Limits	Preview	//Index	History	С	lipboard	Deta	ils
About Entre	<u>z</u>	Display Abstra	ct	Show	y. 20 🗶 So	rt 🤶	Send to	o Text	M

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities


PubMed Services
Journals Database
MeSH Database
Single Citation Matcher
Batch Citation Matcher
Clinical Queries
LinkOut
Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ 1: Proc Inst Mech Eng [H]. 2001;215(6):543-8.

Related Articles, Links



## Contact mechanics of a novel metal-on-metal total hip replacement.

Besong AA, Lee R, Farrar R, Jin ZM.

Department of Mechanical and Medical Engineering, University of Bradford, UK.

The contact mechanics of a novel metal-on-metal total hip replacement (THR) were investigated in this study. The metal-on-metal prosthesis considered consists of a cobalt-chrome acetabular insert connected to a titanium shell through a taper contact, articulating against a cobalt-chrome femoral head. Both the experimental measurement of the displacement of the acetabular insert and the contact area between the two bearing surfaces, and the corresponding numerical predictions using the finite element method have been conducted. Excellent agreement has been demonstrated between the experimental measurement and the finite element prediction under various loads up to 3 kN. The maximum contact pressure at the articulating surfaces has been predicted to be about 31 MPa from a simple axisymmetric finite element model, significantly lower than that of a similar cup but with a monoblock construct. This has been mainly attributed to the flexibility of the insert, leading to an increase in the conformity between the femoral head and the acetabular insert. In addition, the predicted maximum contact pressure is only slightly increased to 37 MPa, from a more realistic three-dimensional anatomical finite element model. The design features on metal-onmetal THRs have been shown to reduce contact stresses and may improve tribological performances of these hard-on-hard bearing couples.

PMID: 11848386 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 ▼ Sort ▼ Send to Text ▼

Write to the Help Desk


NCBI | NLM | NIH

Department of Health & Human Services

Freedom of Information Act | Disclaimer

Apr 19 2004 06:53:45

ALLEL ARRONS -- ... -- 110/0386 5/4/7004







				•					
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	OMIMO	PMC	Journals	Book
Te	PubMed	a for					Go C	lear	
Search L		Limits	Previ	ew/index	History	C	lipboard	Deta	ils
About Entrez		Display Abstra	ct	Sho	w: 20 💆 So	n 🙀	Send t	Text	M

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ **1:** Clin Orthop. 1996 Aug;(329 Suppl):S244-55.

Related Articles, Links

# Metal on metal bearings. A practical alternative to metal on polyethylene total joints?

Black J.

IMN Biomaterials, King of Prussia, PA, USA.

Metal on metal articulation is proposed as an alternative to metal on polymer in total hip replacement arthroplasty as a technical means of reducing wear debris production and subsequent osteolysis leading to the need for surgical revision. The question of whether metal on metal articulation is a practical alternative to current practice is essentially that of whether it is as safe as, and more effective than, metal on polymer articulations in use for more than 20 years. Unfortunately, the metal on metal articulation introduces additional biologic risks associated with production of increased metallic corrosion and wear products. The clinical longevity and success of metal on polymer articulation in total hip replacements, as embodied in the Charnley type, is such that it may prove humanly impossible to determine that metal on metal articulations are more effective, even if that is objectively the case. Therefore, it is suggested that, consistent with modern technical and ethical standards, it cannot be concluded that metal on metal articulation is a practical alternative to current metal on polymer designs. It is suggested that future improvement in total hip replacement arthroplasty outcome is more likely to be through evolutionary than revolutionary designs.

PMID: 8769338 [PubMed - indexed for MEDLINE]

	•
李哲士为《1800·00·00·00·00·00·00·00·00·00·00·00·00·	
Dienlay Abstract	20 Sort Send to Text
Display Abstract	Show: 17

Write to the Help Desk

NCBI | NLM | NIH

Department of Health & Human Services
Freedom of Information Act | Disclaimer

Apr 19 2004 06:53:45


20 2-11 Abstract 0.11ct 111/2-0760338 5/5/2004

PubMed Central

Privacy Policy

Display

Abstract




Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer

Apr 19 2004 06:53:45

5/5/2004

-01:- ...A--0E270E0



Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents **NLM Gateway** TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ 1: J Bone Joint Surg Br. 2002 Jan;84(1):128-36. Adverse tissue reactions to wear particles from Co-alloy

Related Articles, Links

articulations, increased by alumina-blasting particle contamination from cementless Ti-based total hip implants. A report of seven revisions with early failure.

Bohler M, Kanz F, Schwarz B, Steffan I, Walter A, Plenk H Jr, Knahr K.

Institute for Histology and Embryology, Vienna, Austria.


We revised seven alumina-blasted cementless hip prostheses (Ti-alloy stems, cp Ti threaded sockets) with low- or high-carbon Co-alloy bearings at a mean of 20.1 months after implantation because of pain and loosening. Histological examination of the retrieved periprosthetic tissues from two cases in which the implant was stable and three in which the socket was loose showed macrophages with basophilic granules containing metal and alumina wear particles and lymph-cell infiltrates. In one of the two cases of stem loosening the thickened neocapsule also contained definite lymphatic follicles and gross lymphocyte/plasma-cell infiltrates. Spectrometric determination of the concentration of elements in periprosthetic tissues from six cases was compared with that of joint capsules from five control patients undergoing primary hip surgery. In the revisions the mean concentration of implant-relevant elements was 693.85 microg/g dry tissue. In addition to Cr (15.2%), Co (4.3%), and Ti (10.3%), Al was predominant (68.1%) and all concentrations were significantly higher (p < 0.001) than those in the control tissues. The annual rates of linear wear were calculated for six implants. The mean value was 11.1 microm (heads 6.25 microm, inserts 4.82 microm). SEM/EDXA showed numerous fine scratches and deep furrows containing alumina particles in loosened sockets, and stems showed contamination with adhering or impacted alumina particles of between 2 and 50 microm in size.

PMID: 11837818 [PubMed - indexed for MEDLINE]

Dienlay Abstract	20 Sort Send to Text	~
Display Abstract	Show: I	

Write to the Help Desk NCBI NLM NIH Department of Health & Human Services Freedom of Information Act | Disclaimer

Privacy Policy



Display

Abstract


Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer

Sort

Apr 19 2004 06:53:45

Send to

268



Show: 20

Display

Abstract

Sort

Send to

Text

Clinical Queries LinkOut

Related Resources

Order-Documents NLM Gateway

Consumer Health

PubMed Central

Clinical Alerts ClinicalTrials.gov

Privacy Policy

Cubby

TOXNET







S NCB	9	rud	ME	eu	of Me	dicine <b>N</b>	ILM	
Entrez PubMed	Nucleotide	Protein	Genome	Structure	OMIM	PMC	Journals	Books
Search PubMed	for	**************************************				Go C	lear	
	Limits	Previ	ew/Index	History	С	lipboard	De	tails
About Entrez	<b>Display</b> Abstr	ract	Sho	w: 20 👺 Sor	, A	Send to	Text	Ž
Text Version								
Entrez PubMed Overview	🗆 1: J Bone Joi	nt Surg Br.	1997 Mar;79	9(2):316-21.			Related A	Articles, Links
Help   FAQ Tutorial New/Noteworthy E-Utilities	Erratum ir • JB		Gurg Br 1997	Sep;79(5):885.				
PubMed Services Journals Database MeSH Database	Elevate surface		cobalt w	rith metal-o	n-met	al artic	ulating	
Single Citation Matcher	Brodner	W. Bitzan	P. Meisina	er V, Kaider A	, Gottsa	uner-Wo	If F, Kotz i	R.

Department of Occupational Medicine, University of Vienna, Austria.

We determined serum cobalt levels in 55 patients by atomic absorption spectrophotometry before and after implantation of uncemented total hip arthroplasties. In a randomised, prospective trial 27 wrought Co-28Cr-6Mo-0.2C metal-on-metal articulations were compared with 28 ceramic-on-polyethylene hips which did not contain cobalt. Other sources of iatrogenic cobalt loading were excluded. The metal-on-metal group produced detectable serum cobalt levels (median 1.1 microg/l after one year) which were significantly different (p < 0.0001) from those of the ceramic-on-polyethylene control group (median below detection limit of 0.3 microg/l after one year). Our findings indicate that metal-on-metal bearings generate some systemic release of cobalt.

#### Publication Types:

- Clinical Trial
- Randomized Controlled Trial

PMID: 9119865 [PubMed - indexed for MEDLINE]

	20 El Sort	Send to Text
Display Abstract	Show: 20 Son	M Send of Lext M

Write to the Help Desk

NCBI | NLM | NIH


Department of Health & Human Services

Freedom of Information Act | Disclaimer

Jun 1 2004 06:47:16

270

Privacy Policy




the same prosthesis type and no known renal disease. CONCLUSION: Chronic renal failure seems to be responsible for the marked elevation of serum cobalt and serum chromium. CLINICAL RELEVANCE: Despite evidence of adverse health reactions, a possible effect of long-term cobalt and chromium loading cannot be neglected. In our opinion, metal-on-metal bearings in THA should not be inserted in patients with chronic renal failure. Follow-up investigations (serum cobalt, serum chromium, serum creatinine, BUN, echocardiography) should be performed at short intervals.

PMID: 11084743 [PubMed - indexed for MEDLINE]

Show: 20 Send to Text Sort Display Abstract

> Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer



PMID: 10656320 [PubMed - indexed for MEDLINE]

Display Abstract Show: 20 Sort Send to Text


Write to the Help Desk

NCBI | NLM | NIH

Department of Health & Human Services

Freedom of Information Act | Disclaimer

Apr 19 2004 06:53:45







	_	Y		OI 1-TC-MICHIA	
Entrez PubMed	Nucleotide	Protein Ger	nome Structure	OMIM PMC	Journals Books
BubMed	for			Go	Clear
Search Fubined	Limits	Preview/Ind	ex History	Clipboard	Details
About Entrez	Display Abstra	act {	Show. 20	Sort 🔀 Ser	d to ∃Text

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

**PubMed Services** Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ 1: J Arthroplasty. 1991;6 Suppl:S5-10.

Related Articles, Links

### Survivorship analysis of the Ring hip arthroplasty.

Bryant MJ, Mollan RA, Nixon JR.

Department of Orthopaedic Surgery, Musgrave Park Hospital, Belfast, Northern Ireland.

Two hundred fifty-three Ring mark 2 metal-on-metal hip arthroplasties performed between 1968 and 1974 were evaluated using survivorship analysis. Using revision as the criterion for failure, the authors found a cumulative survival rate of 60.4% after 21 years. The results are compared with data from previous studies that used survivorship analysis for metal-on-metal hip arthroplasties, and it is shown that the Ring hip arthroplasty performed as well as the McKee-Farrar prosthesis and better than the Stanmore prosthesis.

PMID: 1774571 [PubMed - indexed for MEDLINE]

Display Abstract	Show:	20 🚼	Sort 🔄	Send to Text	
		**************************************			

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer







				٦.			2110	Journals	Book
<i></i>	PubMed	Nucleotide	Protein	Genome	Structure	OMIM	PMC T generalische granden	PRODUTES	0000
Entrez I							Go Cl	ear	
Search l	PubMed	for L	- \$1497 <del></del>		11:-40		lipboard	Detai	ls
		Limits	Prev	iew/Index	History	Ŭ	11000.0		
bout Entri	ez 🦎					8.80	Send to	Tavi	¥
		Display Abstra	ect	<b>.</b>	Show: 20 😽 Sort		Sendid	rext	20X.01
		# 1000 H 5 10 1 25 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1		THE RESERVE OF THE PARTY OF THE	Dil Dil I Manual Manual Annual Control of the Contr				

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

L 1: Acta Orthop Scand. 2002 Oct;73(5):506-12.

Related Articles, Links

Positive cytokine production in failed metal-on-metal total hip replacements.

Campbell PA, Wang M, Amstutz HC, Goodman SB.

Joint Replacement Institute, Orthopaedic Hospital, Los Angeles, CA 9000, USA. pcampbell@laoh.ucla.edu

Tissues surrounding failed conventional total hips have been shown to produce inflammatory cytokines that can induce osteoclastic bone resorption. We evaluated the cytokine profiles of tissues from 5 failed metal-on-metal total hip replacements. Serial frozen sections were stained using immunohistochemical and in situ hybridization techniques. Inflammatory and osteoclast-stimulating cytokines were noted in the tissues. As compared to a group of 5 metal-polyethylene hip tissues, we found fewer CD68 positive macrophages, and lower levels of TGF-beta and TNF-alpha, but no differences in CD3 positive lymphocytes, IL-1beta, IL-6 and PDGF-alpha in the metal-on-metal tissues. This may be due, in part to the presence of wear particles from sources other than the bearing surfaces. Thus, cytokines associated with bone resorption and implant loosening may occur in total hips despite the use of alternative bearing materials.


PMID: 12440492 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort Send to Text

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer







9.3				•					
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	MIMO	PMC	Journals	Book
Г	PubMed	for			Harris of the same		Go C	lear	
Search   Publica		Limits	Previe	ew/Index	History	CI	ipboard	Detai	ls
About Entre	z	<b>Display</b> Abstra	nct	Sho	w: 20 🖁 Sc	nt 🖁	Send to	Text	

About Entrez

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher **Batch Citation Matcher** Clinical Queries LinkOut Cubby

Related Resources Order Documents **NLM** Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

1: Clin Orthop. 1996 Dec;(333):96-107.

Related Articles, Links

#### Engineering issues and wear performance of metal on metal hip implants.

Chan FW, Bobyn JD, Medley JB, Krygier JJ, Yue S, Tanzer M.

Jo Miller Orthopaedic Research Laboratory, Montreal General Hospital and McGill University, Quebec, Canada.

A major concern in total hip arthroplasty is the generation of polyethylene wear particles at the articulating surfaces and resulting macrophage mediated periimplant osteolysis. There is renewed interest in metal on metal bearings as a solution to this problem in view of their potential for greatly improved wear performance. Using a commercially available hip simulator, the wear performance of metal on metal femoral head and acetabular cup combinations was evaluated and various parameters affecting metal on metal implant wear were identified. Nine implants custom manufactured from 2 medical grades of CoCrMo alloy (ASTM F1537-95 and F75-92) were tested within bovine serum as the lubricant to 3 million cycles (equivalent to approximately 3 years of service in vivo). The progressive wear of the components was determined by gravimetric methods at approximately every 300,000 cycles. The wear rates were characterized by an initial period of accelerated wear after which a lower steady state wear rate was observed for subsequent cycles. The presence of calcium phosphate films on the component surfaces, the microstructure of the lower carbon, wrought alloy, and increased effective radii (decreased diametral clearances) were identified as factors that may be favorable to improved wear performance. The extent of the effect on wear of each parameter, however, cannot be discerned at this point and necessitates a study in which parametric changes are more tightly controlled. The present study suggests that the use of metal on metal articulating surfaces may mitigate the problem of osteolysis by offering improved wear performance.

PMID: 8981885 [PubMed - indexed for MEDLINE]

	<u> </u>	
Display Abstract	Show: 20 Sort Send to Text	¥.

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer







			7	N .		
Entrez	PubMed	Nucleotide	Protein Genom	e Structure	OMIM PMC	Journals Books
Search F	ubMed	for			Go go	lear
Search L		Limits	Preview/Index	History	Clipboard	Details
About Entrez		<b>Display</b> Abstrac	et 👹	Show: 20 Sc	ort ⊱ Sendit	o Text

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy F-I tillities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ **1:** J Bone Joint Surg Br. 2003 Aug;85(6):913-7.

Related Articles, Links

# Levels of metal ions after small- and large-diameter metal-on-metal hip arthroplasty.

Clarke MT, Lee PT, Arora A, Villar RN.

Orthopaedic Surgery, SUNY-Upstate Medical University, Syracuse, New York 13202, USA.

Metal-on-metal (MOM) bearings for hip arthroplasty are increasing in popularity. Concern remains, however, regarding the potential toxicological effects of the metal ions which these bearings release. The serum levels of cobalt and chromium in 22 patients who had undergone MOM resurfacing arthroplasty were compared with a matched group of 22 patients who had undergone 28 mm MOM total hip arthroplasty (THA). At a median of 16 months (7 to 56) after resurfacing arthroplasty, we found the median serum levels of cobalt and chromium to be 38 nmol/l (14 to 44) and 53 nmol/l (23 to 165) respectively. These were significantly greater than the levels after 28 mm MOM THA which were 22 nmol/l (15 to 87, p = 0.021) and 19 nmol/l (2 to 58, p < 0.001) respectively. Since the upper limit for normal patients without implants is typically 5 nmol/l, both groups had significantly raised levels of metal ions. MOM bearings of large diameter, however, result in a greater systemic exposure of cobalt and chromium ions than bearings of small diameter. This may be of relevance for potential long-term side-effects. It is not known to what extent this difference is due to corrosion of the surfaces of the component or of the wear particles produced.

PMID: 12931818 [PubMed - indexed for MEDLINE]

	•	
Display Abstract	Show: 20 Sort Send to Text	3

Write to the Help Desk

NCBI | NLM | NIH

Department of Health & Human Services

Freedom of Information Act | Disclaimer

Help | FAQ Tutorial

New/Noteworthy E-Utilities

PubMed Services Journals Database

Batch Citation Matcher Clinical Queries

Related Resources Order Documents

**NLM Gateway** 

Clinical Alerts

Privacy Policy


ClinicalTrials.gov PubMed Central

TOXNET Consumer Health

MeSH Database Single Citation Matcher

LinkOut

Cubby



Dislocation after total hip replacement in relation to metal-onmetal bearing surfaces.

Clarke MT, Lee PT, Villar RN.

BUPA Cambridge Lea Hospital, Impington, England, UK.

Metal-on-metal (MOM) is a commonly used bearing notable for its 'suction fit' when lubricated. In this study, we examined the capacity for MOM bearings to protect against dislocation after total hip replacement (THR). We undertook a clinical investigation to compare the rate of dislocation of MOM bearings with those of ceramic-on-polyethylene (COP) bearings and found that one MOM bearing dislocated in a series of 109 hips (0.9%) compared with nine of 145 hips (6.2%) in the COP group (p = 0.02). We also performed an in vitro investigation comparing the peak forces generated during forced separation of the two bearings of the same dimensions at velocities from 1 to 50 cm/s. This revealed that the MOM bearing generated significant resistance to separation at all velocities (maximum mean 24 N), whereas the COP did not (maximum mean 1.9 N, p < 0.001). We conclude that MOM bearings are more stable to dislocation than COP bearings as a result of the interfacial forces provided by a thin, lubricating fluid.

PMID: 12892184 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort Send to Text

Write to the Help Desk

NCBI | NLM | NIH

Department of Health & Human Services

Freedom of Information Act | Disclaimer

S NCBI		Pub C	Aed		National Library Medicine	Journals	Books
Entrez PubMed	Nucleotide	T TOTALL			Go Ci	ear	
Search Publied	Limits	Preview/\nd	ex His	itory	Clipboard	Details	
About Entrez	<b>Display</b> Abst		Show: 20	Sort	Send to	Text	
Text Version  Entrez PubMed Overview Help   FAQ Tutorial New/Noteworthy E-Utilities		ab Sci. 1996 Mar risk from or			s.	Related Article	as, Links
PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby  Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central  Privacy Policy	Permaner over the I prosthese Prolonged the most reported; population Combining examined replacem. There was with a total deficits of but combined to but combined the surgicases. All alarming prosthesis malignant large contributions and the type each patal ion time since using ortal term care.	nt replacement of last 50 years. Viges is required. So discontact of metals and they provide direct metals. They provide direct metals shad cancer risk in 40 ent. Overall, the last an 8 percent estal of 347 cases of cancers of the loined results from sk in the first 10 hout there was a gery, with a related three studies who have complicating and composition and composition in cases and composition would be insertion would be insertion would be insertion would be insertion would be insertion.	f joints damaged plance over posture of the metal alloys with better cohost studies. The designation of the prosture of time they like in excess of haemed better well-designation of the prosture of the pros	ed by fractures ible long-tal componer ody fluids redies of persove observation of risk followed up cancer was opoietic malable. 08, 95 persove observation of risk of the cancer was opoietic malable. 08, 95 persove observation of the cancer cancer cancer cases on the cancer cases of th	re or arthritis had been adverse efforts are potential esults in gradual ons with a hip property of these estimates. College for a mean 8.5 and 1.02 (95 perceignancy (leukastent CI 0.97 to re seen in the transitional significant CI 1.00 to 1.1 cuted. Their rest the number of prosthesis will in help quantify and is needed, inclusionated of tissue for and of tissue for could be done of the could be done of the could be done of the done of the could be done of the could	es become corfects of metal lly carcinogeni corrosion of erosthesis have sk in a human studies are sinectively, the st years after him CI 1.00 to 1 mia and lympl 1.20). Signification of the control of	mmon ic. even been milar. cudies p 1.05). homa), ant idies, iity. e years ,005 a gister of d, a ion on res for nt , and ly,

Send to Text

Display Abstract

Show: 20 Sort





National Library of Medicine

	MODI	-	STATE AND PERSONS	A-4		or Me	Of Medicine		
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	MIMO	PMC	Journals	Book
Search F	<sup>2</sup> ubMed	for					Go CI	ear	
Scarcition	**************************************	Limits	Previ	ew/Index	History	C	lipboard		etails
About Entrez		<b>Display</b> Abstra	ct	Show	w: 20 😮 Sc	ort 🕺	Send to	Text	

Text Version

Entrez PubMed Overview Heip | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services
Journals Database
MeSH Database
Single Citation Matcher
Batch Citation Matcher
Clinical Queries
LinkOut
Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ **1:** J Bone Joint Surg Br. 2004 Mar;86(2):177-84.

Related Articles, Links

# Metal-on-metal resurfacing of the hip in patients under the age of 55 years with osteoarthritis.

Daniel J, Pynsent PB, McMinn DJ.

The Birmingham Nuffield Hospital, Edgbaston, Birmingham, England, UK.

The results of conventional hip replacement in young patients with osteoarthritis have not been encouraging even with improvements in the techniques of fixation and in the bearing surfaces. Modern metal-on-metal hip resurfacing was introduced as a less invasive method of joint reconstruction for this particular group. This is a series of 446 hip resurfacings (384 patients) performed by one of the authors (DJWM) using cemented femoral components and hydroxyapatite-coated uncemented acetabular components with a maximum follow-up of 8.2 years (mean 3.3). Their survival rate, Oxford hip scores and activity levels are reviewed. Six patients died due to unrelated causes. There was one revision (0.02%) out of 440 hips. The mean Oxford score of the surviving 439 hips is 13.5. None of the patients were told to change their activities at work or leisure; 31% of the men with unilateral resurfacings and 28% with bilateral resurfacings were involved in jobs that they considered heavy or moderately heavy; 92% of men with unilateral hip resurfacings and 87% of the whole group participate in leisure-time sporting activity. The extremely low rate of failure in spite of the resumption of high level occupational and leisure activities provides early evidence of the suitability of this procedure for young and active patients with arthritis.

PMID: 15046429 [PubMed - indexed for MEDLINE]


A CONTRACTOR AND A CONT	2.00 PM
	- 155 2011 C : 3 (2) 20 20 20 20 31 T A 24 (2) 1 (2)
Display Abstract	120 \$1 Son \$2 8 Sending 1 ext \$1
Citylopidy Mostract	Show: 20
	SHUW.

Write to the Help Desk

NCBI | NLM | NIH

Department of Health & Human Services

Freedom of Information Act | Disclaimer



Write to the Help Desk

NCBI | NLM | NIH

Department of Health & Human Services

Freedom of Information Act | Disclaimer

Sort

Show: 20

Display

Abstract

Jun 1 2004 06:47:16

Send to

Text







<u> </u>			*	•			
Entrez	PubMed	Nucleotide	Protein Gend	ome Structure	OMIM P	PMC Journals Boo	١
Search P	ubMed	for			G	o Clear	
Journal III		Limits	Preview/Inde	x History	y Clipbo	eard Details	
About Entrez		<b>Display</b> Abstra	ct E	Show: 20	Sort 😿 S	end to Text	

About Entrez

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents **NLM Gateway** TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

**1:** Clin Orthop. 1996 Aug;(329 Suppl):S206-16.

Related Articles, Links

#### Metal versus polyethylene wear particles in total hip replacements. A review.

Doorn PF, Campbell PA, Amstutz HC.

Joint Replacement Institute at Orthopaedic Hospital, Los Angeles, CA, USA.

Research has recently been focused on the development of hip replacements with alternative bearing surfaces with cobalt chrome alloy, to avoid the production of polyethylene wear particles in hip replacements and polyethylene wear debris mediated bone lysis. Cobalt chrome on cobalt chrome bearing surfaces are being reevaluated. Characterization of wear particles and studies on the reaction of the body to these particles, have played an important role in the determination of the factors that cause aseptic loosening and will therefore play an important role in the comparison of metal on polyethylene and metal on metal hip prostheses. In this paper, a comparison between the different aspects of metal and polyethylene wear particles is made using data from the literature and the authors' experience. The authors conclude that techniques need to be optimized to isolate and characterize individual metal wear particles from periprosthetic tissues and they advocate the performance of in vitro studies with these in vivo generated wear particles or comparable particles.


#### Publication Types:

- Review
- Review, Tutorial


PMID: 8769335 [PubMed - indexed for MEDLINE]

Display Abstract	Chaum	20	$\overline{\mathbf{Y}}$	Sort	<b>.</b>	Send to	Text	×
A Section of the Control of the Cont	Show:		البنت					

Write to the Help Desk NCBI NLM NIH Department of Health & Human Services Freedom of Information Act | Disclaimer







			Y		04 7/254	M17-1120	<del></del>	
Entrez	PubMed	Nucleotide	Protein Geno	ome Structure	ОМІМ	PMC	Journals	Book
Search	PubMed	for				Go Cl	ear	
ocarca,	Assessment of the state of the	Limits	Preview/Inde	x History	CI	ipboard	Deta	ils
About Entre	e <b>z</b>	<b>Display</b> Abstrac	ot E	Show: 20	Sort 😸	Send to	Text	

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents **NLM Gateway** TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ 1: J Biomed Mater Res. 1998 Oct;42(1):103-11. **SWILLTY** 

InterScience

Related Articles, Links

Metal wear particle characterization from metal on metal total hip replacements: transmission electron microscopy study of periprosthetic tissues and isolated particles.

Doorn PF, Campbell PA, Worrall J, Benya PD, McKellop HA, Amstutz HC.

Joint Replacement Institute at Orthopaedic Hospital, Los Angeles, California 90007, USA.

The less intense tissue reaction around metal on metal total hip replacements (THRs) compared to metal on polyethylene (PE) THRs may be explained by the differences in the characteristics of metal wear particles. In this study, transmission electron microscopy was used to study metal wear particles that were either in situ in cells or had been extracted from the cells by a new technique based on enzymatic tissue digestion. The tissues were obtained from 13 patients undergoing revision of metal on metal THRs with cobalt-chromium-molybdenum (CoCrMo) bearing couples. Most of the CoCrMo wear particles were smaller than 50 nm (range 6-834 nm) and round to oval in shape with irregular boundaries. This size range is considerably smaller than that reported for PE particles. While even a small volume of metal wear will produce high numbers of particles, the apparently less severe local tissue reaction to metal particles may be due to the possibility that corrosion, dissolution, and dissemination of metal particles may result in fewer local biological effects than the long-term retention of PE particles in the periprosthetic tissues.

PMID: 9740012 [PubMed - indexed for MEDLINE]

Display Abstract	4 T. 150 (\$2120)[	3
Display Abstract X		<u> </u>

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer







Entrez	PubMed	Nucleotide F	Protein Genor	re Structure	OMIM PMC	Journals Book
Search Pu	bMed	for	······································	Harris Marie Composition of the	(Go	Clear
		Limits	Preview/Index	History	Clipboard	d Details
About Entrez		Display Abstrac	t I	Show: 20 🔄	Sort 💆 Ser	nd to Text

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ **1:** Clin Orthop. 1996 Aug;(329 Suppl):S187-205.

Related Articles, Links

#### Tissue reaction to metal on metal total hip prostheses.

Doorn PF, Mirra JM, Campbell PA, Amstutz HC.

Joint Replacement Institute, Los Angeles, CA, USA.

The periprosthetic tissue reaction to polyethylene wear debris in metal on polyethylene total hip replacements is strongly implicated as the cause of osteolysis. This has led to a renewed interest in metal on metal total hip replacements. However, little is known about the role of wear debris in failures of these prostheses. Capsular and interface tissues from 9 long and short term metal on metal total hip replacement retrievals were studied to assess the tissue reaction around these prostheses. As compared with metal on polyethylene cases, the extent of the granulomatous inflammatory reaction and the presence of foreign body type giant cells was much less intense in metal on metal cases, likely because of the lower numbers and overall smaller size of metal wear debris particles. This may lead to a better transport of the particles from the joint tissues and a lower incidence of periprosthetic osteolysis around metal on metal hip replacement.

PMID: 8769334 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort Send to Text
------------------	----------------------------

Write to the Help Desk

NOB | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer

Entrez PubMed Overview Help | FAQ

PubMed Services

Journals Database

MeSH Database Single Citation Matcher Batch Citation Matcher

Clinical Queries

Related Resources

Order Documents

Consumer Health

PubMed Central

NLM Gateway TOXNET

Clinical Alerts Clinical Trials.gov


Privacy Policy

LinkOut Cubby

Tutorial New/Noteworthy E-Utilities







				•	01 1 100101110 -	
Entrez	PubMed	Nucleotide	Protein Genor	ne Structure	OMIM PMC	Journals Books
Search F	PubMed	for		and the second s	Go	lear
Journal I	on training the same	Limits	Preview/Index	History	Clipboard	Details
About Entrea	?	<b>Display</b> Abstra	ct 💆	Show: 20 So	ort 🕌 Sendit	Text
Text Version	ı	<del>-</del>				

☐ **1:** Clin Orthop. 1996 Dec;(333):108-17.

Related Articles, Links

#### Modern metal on metal articulation for total hip replacements.

Dorr LD, Hilton KR, Wan Z, Markovich GD, Bloebaum R.

Center for Arthritis and Joint Implant Surgery, University of Southern California University Hospital, Los Angeles, CA 90033, USA.

Between 1991 and 1994, 70 patients received total hip replacements with metal on metal articulation. The results of 54 of these patients with 54 hips who have a 2- to 4year (2.7-year average) followup are reported. Patients were prospectively evaluated using the Harris hip score, a patient self assessment form, and radiographs. Hip aspiration was performed preoperatively and 6 to 24 months postoperatively in 24 hips with metal on metal articulations. Implant retrieval was obtained from 2 patients. Harris hip score averages increased from 49 to 93. No patient had revision surgery for loosening, but 1 had revision surgery for dislocation. Patient self assessment forms showed 51 of 54 patients scored their results as good or excellent. Serial radiographs did not show loosening or osteolysis. Wear could not be measured radiographically. Synovial fluid samples had metal particles of 1 to 10 microm in 10 hips. Twenty patients had bilateral total hip replacements with 1 hip metal on polyethylene articulation, and patients could not determine any difference between the hips. Compared with historic results of previous metal on metal prostheses, the modern metal on metal articulation investigated in this study did not have early acetabular loosening or clinical symptoms of component impaction. Retrieval implants and synovial fluid analysis suggest early wear was minimal.

PMID: 8981886 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort Send to Text


Write to the Help Desk

NCBI | NLM | NIH

Department of Health & Human Services
Freedom of Information Act | Disclaimer







······································									
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	MIMO	PMC	Journals	Books
[F	ubMed	for					Go C	lear	
Search L		Limits	Previ	ew/Index	History	C	lipboard	Deta	ils
About Entrez		Display Abstra	ct	Sho	w: 20 🔻 So	rt 😽	Send to	Text	Z

Entrez PubMed Overview Heip | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ **1:** J Bone Joint Surg Am. 2000 Jun;82(6):789-98.

Related Articles, Links

Comment in:

• J Bone Joint Surg Am. 2001 May;83-A(5):783-5.

Full text article at www.cjbjs.org

Total hip arthroplasty with use of the Metasul metal-on-metal articulation. Four to seven-year results.

Dorr LD, Wan Z, Longjohn DB, Dubois B, Murken R.

Weber Institute, St. Gallen, Switzerland.


BACKGROUND: Total hip replacements with a metal-on-metal articulation were commonly used until the mid-1970s; most were then abandoned in favor of hip replacement with a metal-on-polyethylene articulation. The reason for this change was primarily early cup loosening, which was more prevalent with these metal-on-metal designs than it was with metal-on-polyethylene designs. In the late 1980s, a metal-onmetal design with improved clearance (adequate space between the femoral head and the acetabular articulation surface to allow fluid film lubrication and clearance of any debris from within this joint), metal hardness, and reproducible surfaces was introduced by Sulzer Orthopedics in Switzerland. Orthopaedic surgeons were interested in this Metasul articulation because the contribution of polyethylene wear particles to the failure of total hip replacements had become evident. This study was undertaken to review the clinical performance of this implant and to determine if early acetabular loosening or revision and wear and osteolysis were prevalent. METHODS: Between 1991 and 1994, seventy patients (seventy hips) had a total hip replacement with the Metasul metal-on-metal articulation and a cemented Weber cup. Nine patients died less than four years after the replacement; none of these deaths were related to the operation. Five patients were not available for radiographic evaluation, but they were contacted and it was known that the hip was not painful and had not been revised. Fifty-six patients (fifty-six hips) had complete clinical and radiographic data four to 6.8 years after the operation, and they made up the study group. The patients were evaluated with use of the Harris hip score, a patient-self-assessment form, and radiographs. RESULTS: At an average of 5.2 years (range, four to 6.8 years) after the operation, the average total Harris hip score for the fifty-three patients who did not have a revision was 89.6 points (range, 62 to 100 points). The average Harris pain score was 41.0 points (range, 30 to 44 points), and the average Harris limp score was 9.4 points (range, 5 to 11 points). One patient had revision of a loose cup, but there were no other loose acetabular components in the series. Two patients had revision of the acetabular component because of dislocation. No patient had a loose or revised femoral component. Therefore, the mechanical failure rate was one (2 percent) of fifty-six patients. Thirty-six of forty-seven patients who completed the patient-self-assessment form rated their result as excellent; seven, as very good; two, as good; one, as fair; and one, as poor. Wear could not be measured on radiographs because of the metal-onmetal articulation. No hip had radiographic evidence of acetabular osteolysis and two

hips had calcar resorption, but there was no other radiographic evidence of focal osteolysis. CONCLUSIONS: Our four to seven-year experience with this articulation surface indicates that the clinical results are similar to those of total hip replacements with a metal-on-polyethylene articulation. We believe that the Metasul articulation may have a role in reducing the wear that occurs with total hip replacement. The Metasul articulation appears to be particularly indicated for more active patients. A historical comparison with the reports in the literature of which we are aware indicated that the hips in our study had a lower rate of acetabular revision and loosening than did those with previous metal-on-metal designs and that they had no more acetabular loosening or osteolysis than did those with metal-on-polyethylene articulations followed for an average of five years.

PMID: 10859098 [PubMed - indexed for MEDLINE]

Annual desirability of the control o	Sort Send to Text	
Dienlay Abstract		
	Show:	

Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer



Related Resources
Order Documents
NLM Gateway
TOXNET
Consumer Health
Clinical Alerts
ClinicalTrials.gov
PubMed Central


Privacy Policy

Write to the Help Desk

NCBI | NLM | NIH

Department of Health & Human Services
Freedom of Information Act | Disclaimer

Privacy Policy




PMID: 3403594 [PubMed - indexed for MEDLINE]

Display Abstract Show: 20 - Sort Send to Text

Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer







				•					
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	MIMO	PMC	Journals	Books
Search F	oubMed	for					Go C	lear	
Search	THE RESERVE OF THE PARTY OF THE	Limits	Previ	ew/index	History	С	lipboard	De	etails

About Entrez

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

1: Rev Chir Orthop Reparatrice Appar Mot. 1999 Jun;85(3):217-25.

Show: 20

•85(3)•217-25 Related Articles. Links

Text

Send to

full text on E2MED

Subscription required

## [Value and limits of determining serum cobalt levels in patients with metal on metal articulating prostheses]

Sort

[Article in French]

Abstract

Display

Gleizes V, Poupon J, Lazennec JY, Chamberlin B, Saillant G.

Service de Chirurgie Orthopedique, Traumatologique et Reparatrice de l'Appareil Locomoteur, CHU Pitie-Salpetriere, Universite Paris VI.

PURPOSE OF THE STUDY: The purpose of this study was to measure the serum cobalt levels and their correlation with clinical and radiological findings in patients with metal on metal hip articulating surfaces. METHOD: Forty-one patients with metal on metal hip arthroplasty were reviewed retrospectively at mean follow-up of 12.9 months. Serum cobalt levels were determined for each patient by atomic absorption spectrometry at the maximal follow-up and were compared to a control group (19 patients). Two patients and one control subjects also performed exercise on a treadmill in order to appreciate the influence of physical activity on serum cobalt levels. RESULTS: The metal on metal group presented higher serum cobalt levels than those of the control group (p < 0.0001). There was no correlation between serum cobalt and clinical and radiological findings at the exception of patient age (n = 40, r = 0.37). However, when the followup was greater than 18 months, mean serum cobalt was significantly higher compared to a follow-up less than 18 months. The physical exercise test led to a moderate elevation (around 10 p. 100) of cobalt in the two patients but not in the control subject. DISCUSSION AND CONCLUSION: The interpretation of an elevated cobalt serum level is difficult. Cobalt-containing drugs, other implants, excess of activity and diseases (renal failure) may influence serum cobalt level. In this study, the high serum cobalt levels seem not linked to a failure of the implant, mainly because of the short follow-up. They could rather be attributed to an increase of the patient activity beginning 18 months after the surgery. Because potential long-term cobalt toxicity and carcinogenicity is not well known, careful medical follow-up should be emphasized specially in young patients.

PMID: 10422126 [PubMed - indexed for MEDLINE]

	il and the second secon	Appart of the Control of
Display Abstract	Su 20 Sort Se	end to Text
Display Abstract	Show:	21.7 (1.9 2. 3223)







	D. Jankani	Nucleotide	Protein Geno	me Structure	MIMO s	PMC	,loumais	Books
Entrez	PubMed					Go Ci	ear	
Search	PubMed	for Limits	Preview/Inde	x Histor	У	Clipboard	Details	s
About Entr	ez	Display Abstrac	et 🖁	Show: 20	Sort	Send to	Text	3

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services
Journals Database
MeSH Database
Single Citation Matcher
Batch Citation Matcher
Clinical Queries
LinkOut
Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

 $\square$  1: Proc Inst Mech Eng [H]. 2000;214(1):39-47.

ingenta select Related Articles, Links

\*\* . \_ \_ ... ... ... ... 40740000 E/N/2004

A comparative joint simulator study of the wear of metal-onmetal and alternative material combinations in hip replacements.

Goldsmith AA, Dowson D, Isaac GH, Lancaster JG.

DePuy International, A Johnson and Johnson Company, Leeds, UK.

While total hip replacement represents the major success story in orthopaedic surgery in the twentieth century, there is much interest in extending even further, early in the twenty first century, the life of implants. Osteolysis has been identified as a major factor limiting the life of prostheses, with indications that fine polyethylene wear debris, generated primarily at the interface between the femoral head and the acetabular cup, promotes the process. There is therefore considerable interest in the introduction of alternative wear resistant systems to limit the deleterious effects of wear. These alternatives include ceramic-on-ceramic and metal-on-metal configurations and the present paper is primarily concerned with the latter. Some six pairs of new metal-onmetal implants of 36 mm diameter and four pairs of existing metal-on-metal implants of 28 mm diameter were tested in a ten-station hip joint simulator in the presence of a 25 per cent bovine serum solution. The implants were tested in the anatomical position to  $5 \times 10(6)$  cycles. The new heads and cups were manufactured from CoCrMo alloy with careful attention being paid to sphericity and surface finish of both components. The wear performance of the new and existing metal-on-metal total hip replacements have been evaluated and compared. The overall wear rates have then been compared with previously reported wear rates for a zirconia-on-polyethylene prosthesis of 22 mm diameter tested on the same simulator. The comparison is taken further by recalling published penetration data for metal-on-polyethylene implants of 22 and 28 mm diameter and converting these to volumetric wear rates. It was found that the heads and cups in metal-on-metal joints wore by almost equal amounts and that the opposing surfaces converged to similar surface roughness as the testing time increased. Steady state wear rates were generally achieved after 1-2  $\times$  10(6) cycles. The mean long-term wear rates for the metal-on-metal prostheses were very low, being 0.36 mm3/10(6) cycles and 0.45 mm3/10(6) cycles for the new implants of 36 mm diameter and established implants of 28 mm diameter respectively. These wear rates compare with 6.3 mm3/10(6) cycles for zirconia-on-ultra-high molecular weight polyethylene tested on the same simulator and representative clinical values for metal-on-polyethylene of 36 mm3/year for heads of 22 mm diameter and a reported range of 60-180 mm3/year for 28 mm heads. These values do not translate directly into numbers of particles, since the metallic debris from metal-on-metal joints is very fine. The number of metallic particles may exceed the number of polyethylene wear particles from an otherwise similar metalon-polyethylene joint by a factor of 10(3). A detailed discussion of the size and morphology of wear debris and tissue reaction to various forms of debris is beyond the scope of this paper, but the biological response to polymeric, metallic and ceramic wear 290 debris forms a major subject for further study. The present investigation nevertheless confirms the potential of carefully designed and manufactured metal-on-metal total replacement joints for the treatment of diseased and damaged hips.

PMID: 10718049 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort Send to Text
	S11(1)(A)

Write to the Help Desk

NCBI | NLM | NIH

Department of Health & Human Services

Freedom of Information Act | Disclaimer

Entrez PubMed

New/Noteworthy

PubMed Services
Journals Database

Batch Citation Matcher Clinical Queries

Related Resources

Order Documents

**NLM Gateway** 

MeSH Database Single Citation Matcher


Overview

Help | FAQ Tutorial


E-Utilities

LinkOut

Cubby







	ستم عدين			7					
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	OMIMO	PMC - #2#####	Journals	Вι
Search	PubMed	for		entagen various de l'accept per especient Navaen demance de desert	epyholegamuun (sin esimpee allan as 1 julgalain maisalesiseen parkiseen kuulkassa.	emilies conjuniementalistic rengellato con	_ Go] <u>⊭C</u>	lear	
,	Services (CCC)	Limits	Previ	ew/Index	History	CI	ipboard	Deta	IIS
		ing \$100 的复数数数数数数数数数数数		Price	00 000	+ 83	Send to	Text	÷
About Ent	rez	Display Abstra	act	Sho	w: 20 💽 Soi	t 🔄	OCI10 N		

Text Version

1: J Bone Joint Surg Am. 2001 Mar;83-A(3):428-36.

Related Articles, Links

Full text article at www.cjbjs.org

Metal sensitivity in patients with orthopaedic implants.

Hallab N, Merritt K, Jacobs JJ.

Department of Orthopaedic Surgery, Rush-Presbyterian-St. Luke's Medical Center, Chicago, IL 60612, USA. nhallab@rush.edu

All metals in contact with biological systems undergo corrosion. This electrochemical process leads to the formation of metal ions, which may activate the immune system by forming complexes with endogenous proteins. Implant degradation products have been shown to be associated with dermatitis, urticaria, and vasculitis. If cutaneous signs of an allergic response appear after implantation of a metal device, metal sensitivity should be considered. Currently, there is no generally accepted test for the clinical determination of metal hypersensitivity to implanted devices. The prevalence of dermal sensitivity in patients with a joint replacement device, particularly those with a failed implant, is substantially higher than that in the general population. Until the roles of delayed hypersensitivity and humoral immune responses to metallic orthopaedic implants are more clearly defined, the risk to patients may be considered minimal. It is currently unclear whether metal sensitivity is a contributing factor to implant failure.

TOXNET
Consumer Health
Clinical Alerts
ClinicalTrials.gov
PubMed Central

Privacy Policy

Publication Types:

Review

Review, Tutorial

PMID: 11263649 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort	Send to Text

Write to the Help Desk

NCBI | NLM | NIH

Department of Health & Human Services

Freedom of Information Act | Disclaimer

Books

8	NCBI
Entrez	PubMed
Search Pi	ıbMed



Preview/Index

. I EN /3

Clipboard

L	ibrary dicine	
MIMO	РМС	Journals
	Go	Clear

Send to

About Entrez

**Text Version** 

Entrez PubMed Overview Help [FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents **NLM Gateway** TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

**1:** J Orthop Res. 2004 Mar;22(2):250-9.

Nucleotide

Display

Limits

Abstract

ELSEVIER

Related Articles, Links

Details

Text

FULL-TEXT ARTICLE Differences in the fretting corrosion of metal-metal and ceramic-metal modular junctions of total hip replacements.

Structure

History

Sort

Hallab NJ, Messina C, Skipor A, Jacobs JJ.

Department of Orthopedic Surgery, Rush Presbyterian St. Lukes Medical Center, 1653 W Congress Parkway, Chicago, IL 60612, USA. nhallab@rush.edu

The use of modular interlocking components is a central design feature of total joint replacements. In this investigation we hypothesized that clinically available ceramicmetal modular connections used in total hip arthroplasty release more metal through fretting corrosion than traditional metal-metal modular connections. This was investigated using an in vitro comparison of ceramic (zirconia, ZrO2) and metal (Coalloy) femoral-head fretting upon Co-alloy stem components. In vitro fretting corrosion testing consisted of potentiodynamic monitoring and analysis of metal release from zirconia and Co-alloy 28 mm femoral heads with similar surface roughnesses (Ra=0.46 microm) on identical Co-alloy stems at 2.2 kN for 1x10(6) cycles at 2 Hz. In contrast to our original hypothesis, we found greater metal release (approximately 11-fold increase in Co and 3-fold increase in Cr) and potentiodynamic fretting of metal-metal modular junctions when compared to ceramic-metal. Potentiodynamic testing demonstrated that lower initial voltages (-266<153 mV), greater maximum voltage changes (116>56 mV, p<0.05, t-test) and voltage variability (3>0.5 mV, p<0.05, t-test) were associated with the open circuit potentials of Co-alloy on Co-alloy junctions when compared to zirconia on Co-alloy junctions. In this study of a single total hip replacement stem and head design, zirconia heads mated with Co-alloy stems produced less fretting than Co-alloy heads mated with Co-alloy stems. Although further studies are necessary with a variety of implant designs and under different experimental conditions, the evidence presented here should, in part, alleviate concerns of increases in fretting corrosion at modular junctions of ceramic-metal coupled components.

PMID: 15013082 [PubMed - indexed for MEDLINE]

And the state of t	ि स्वा ि न	Condito ITov
<b>■Uisplay</b> Abstract	<sub>сь</sub>   20 🔚   50п 🛚 🚎	Send to lext
100	Show: I	

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer







Entrez	PubMed	Nucleotide	Protein	Genome	Structure	OMIM	PMC	Journals	В
Search	PubMed	ubMed for		ing processors were supported to the contract of the contract			Go		
,	Cupillian and Manager or Google of Action for the A	Limits	Previ	ew/Index	History	CI	ipboard	Deta	ails
ль				120	20 🕞 So		Send to	Text	-
About Ent	rez	Display Abstr	act	<b>▼</b> Sho	w: 20 🕝 So	n 🔟	BEIGHT STATE		

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services
Journals Database
MeSH Database
Single Citation Matcher
Batch Citation Matcher
Clinical Queries
LinkOut
Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

□ 1; J Biomed Mater Res. 2000 Sep;53(5):480-9.

Related Articles, Links




A triple assay technique for the evaluation of metal-induced, delayed-type hypersensitivity responses in patients with or receiving total joint arthroplasty.

Hallab NJ, Mikecz K, Jacobs JJ.

Department of Orthopedic Surgery, Rush-Presbyterian-St. Luke's Medical Center, Chicago, IL 60612, USA. nhallab@rush.edu

The determination of biocompatibility has been dominated historically by the characterization of candidate materials based upon the observation of adverse host responses. However, some adverse responses are subtle in clinical settings and continue to foster debate and investigation. One of these responses is "metal allergy" or hypersensitivity to metallic biomaterials. Current methods used to diagnose hypersensitivity reactions, such as dermal patch testing and migration inhibition assays, are not well accepted in orthopedic practice as a means for the characterization of hypersensitivity to metallic joint-replacement components. An increasing need to resolve whether metal sensitivity may be a significant and/or predisposing factor for eliciting an over-aggressive immune response in patients with metallic implant components requires improved and standardized widespread study. Here we present three in vitro methodologies: (1) a proliferation assay, (2) cytokine analysis using ELISA, and (3) a migration inhibition assay. When in conjunction with one another, these assays may be used to more comprehensively quantify metal-induced hypersensitivity responses. Therefore, these methodologies are detailed with the intent of facilitating multi-center large-scale studies. In the following cases, a multi-assay approach for measuring the prevalence of delayed-type hypersensitivity in orthopedic patients shows the propensity to yield a more comprehensive and, therefore, more conclusive determination than currently employed patch testing or single assay techniques. Copyright 2000 John Wiley & Sons, Inc.

PMID: 10984695 [PubMed - indexed for MEDLINE]



Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy

E-Utilities

PubMed Services
Journals Database
MeSH Database
Single Citation Matcher
Batch Citation Matcher
Clinical Queries
LinkOut
Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

Orthopaedic implant related metal toxicity in terms of human lymphocyte reactivity to metal-protein complexes produced from cobalt-base and titanium-base implant alloy degradation.


Hallab NJ, Mikecz K, Vermes C, Skipor A, Jacobs JJ.

Department of Orthopedic Surgery, Rush-Presbyterian-St. Lukes Medical Center, Chicago, IL, USA.

Metal toxicity from sources such as orthopaedic implants was investigated in terms of immune system hyper-reactivity to metal implant alloy degradation products. Lymphocyte response to serum protein complexed with metal from implant alloy degradation was investigated in this in vitro study using primary human lymphocytes from healthy volunteers (n = 10). Cobalt chromium molybdenum alloy (Co-Cr-Mo, ASTM F-75) and titanium alloy (Ti-6Al-4V, ASTM F-136) beads (70 microm) were incubated in agitated human serum at 37 degrees Celsius to simulate naturally occurring metal implant alloy degradation processes. Particulate free serum samples, which were incubated with metal, were then separated into molecular weight based fractions. The amounts of soluble Cr and Ti within each serum fraction were measured and correlated with lymphocyte proliferation response to the individual serum fractions. Lymphocytes from each subject were cultured with 11 autologous molecular weight based serum fractions either with or without added metal. Two molecular weight ranges of human serum proteins were associated with the binding of Cr and Ti from Co-Cr-Mo and Ti implant alloy degradation (at < 30 and 180-330 kDa). High molecular weight serum proteins (approximately 180 kDa) demonstrated greater lymphocyte reactivity when complexed with metal released from Co-Cr-Mo alloy and Ti alloy than with low (5-30 kDa) and midrange (30-77 kDa) serum proteins. When the amount of lymphocyte stimulation was normalized to both the moles of metal and the moles of protein within each fraction (Metal-Protein Complex Reactivity Index, MPCRI), Cr from Co-Cr-Mo alloy degradation demonstrated approximately 10 fold greater reactivity than Ti in the higher molecular weight serum proteins (approximately 180-250 kDa). This in vitro study demonstrated a lymphocyte proliferative response to both Co-Cr-Mo and Ti alloy metalloprotein degradation products. This response was greatest when the metals were complexed with high molecular weight proteins, and with







Entrez	PubMed	Nucleatide	Protein	Genome	Structure	ОМІМ	РМС	Journals	В
Search F	PubMed	for			Conscious demonstration de array opposite autorio de una codo protecido definimientos que se	The same of the sa	Go C	lear	
Format   Format	the state of the s	Limits	Previe	w/Index	History	Cli	ipboard	D	etails
About Entre	:Z	Display Abstra	ct	Shov	w: 20 🔄 Soi	t 🔄	Send to	Text	<u> </u>

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services
Journals Database
MeSH Database
Single Citation Matcher
Batch Citation Matcher
Clinical Queries
LinkOut
Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ 1: J Biomed Mater Res. 2001 Sep 5;56(3):427-36.

InterScience'


Related Articles, Links

Differential lymphocyte reactivity to serum-derived metalprotein complexes produced from cobalt-based and titaniumbased implant alloy degradation.

Hallab NJ, Mikecz K, Vermes C, Skipor A, Jacobs JJ.

Department of Orthopedic Surgery, Rush-Presbyterian-St. Lukes Medical Center, 1653 West Congress Avenue, Chicago, IL 60612, USA. nhallab@rush.edu

The lymphocyte response to serum protein complexed with metal from implant alloy degradation was investigated in this in vitro study using primary human lymphocytes from healthy volunteers (n = 10). Cobalt chromium molybdenum alloy (Co-Cr-Mo, ASTM F-75) and titanium alloy (Ti-6Al-4V, ASTM F-136) beads (70 microm) were incubated in agitated human serum at 37 degrees C to simulate naturally occurring metal implant alloy degradation processes. Particulate free serum samples that had been incubated with metal were then separated into molecular weight based fractions. The amounts of soluble Cr and Ti within each serum fraction were measured and correlated with lymphocyte proliferation response to the individual serum fractions. Lymphocytes from each subject were cultured with 11 autologous molecular weight based serum fractions either with or without added metal. Two molecular weight ranges of human serum proteins were associated with the binding of Cr and Ti from Co-Cr-Mo and Ti implant alloy degradation (at <30 and 180-250 kDa). High molecular weight serum proteins (approximately 180 kDa) demonstrated greater lymphocyte reactivity when complexed with Cr alloy and Ti alloy than low (5-30 kDa) and midrange (30-77 kDa) serum proteins. When the amount of lymphocyte stimulation was normalized to both the moles of metal and the moles of protein within each fraction (metal-protein complex reactivity index), Cr from Co-Cr-Mo alloy degradation demonstrated approximately 10-fold greater reactivity than Ti in the higher molecular weight serum proteins (approximately 180 kDa). This in vitro study demonstrated a lymphocyte proliferative response to both Co-Cr-Mo and Ti alloy metalloprotein degradation products. This response was greatest when the metals were complexed with high molecular weight proteins, and with metal-protein complexes formed from Co-Cr-Mo alloy degradation.







Entrez	PubMed	Nucleotide	Protein	Genome	Structure	OMIM	PMC	Journals	51
Search	PubMed	for Limits	Previ	ew/Index	History	Ci	Go C pboard	N. 12000 N. 1200	tails
About Entr	ez	Display Abstra	act	Show	w: 20 ☑ Sort	Berlinding on a ray 407 (Station	Send to	Text	<u></u>
Text Version	วก	□1: J Biomed	Mater Re	es. 2002 Jun	5;60(3):420-3	3.	Rela	ated Articles,	Links

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services
Journals Database
MeSH Database
Single Citation Matcher
Batch Citation Matcher
Clinical Queries
LinkOut
Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

Concentration- and composition-dependent effects of metal ions on human MG-63 osteoblasts.

Hallab NJ, Vermes C, Messina C, Roebuck KA, Glant TT, Jacobs JJ.

Department of Orthopedic Surgery, Rush-Presbyterian St. Lukes Medical Center, Chicago, Illinois 60612, USA. nhallab@rush.edu

Metal debris from implants has been shown to alter the function of osteoblasts in cell cultures. Its remains unclear, however, if specific forms of released ionic metals are involved in the pathogenesis of periprosthetic osteolysis. We evaluated the relative effects of ionic forms of implant metals by treating human osteoblast-like MG-63 osteosarcoma cells with eight concentrations (0.001-10.0 mM) of Cr(+3), Mo(+5), Al(+3), Ta(+5), Co (+2), Ni(+2), Fe(+3), Cu(+2), Mn(+2), Mg(+2), Na(+2), and V(+3) chloride solutions. The results demonstrated that the metal ions differentially affected osteoblast proliferation, viability, type-I collagen gene expression, and cytokine release. The metal ions were ranked in order from least to most toxic (based on a 50% reduction in viability) as follows: Na < Cr < Mg <Mo < Al < Ta < Co < Ni < Fe < Cu < Mn < V. Metal-induced decreases in osteoblast proliferation were similar in ranking. Nontoxic concentrations of metals had no effect on procollagen alpha1[I] gene expression; only at toxic concentrations did metals produce a decrease in gene expression. The most toxic metals (V, Mn, Fe, and Ni) were also the only metals found to induce IL-6 secretion on a per cell basis (of the cytokines tested, interleukin 6 (IL-6), interleukin beta 1 (IL-1beta), transforming growth factor beta 1 (TGFbeta1), and tumor necrosis factor alpha (TNF-alpha), only IL-6 was detectable in the culture medium after 48 h for any metal at any concentration). Less toxic metals (e.g., Co and Cr) had little effect on IL-6 release, even at high concentrations. In general, metal ions reduced osteoblast function (i.e., proliferation and collagen gene expression) in proportion to the degree of toxicity. These results support the hypothesis that adverse local cellular responses (particularly necrotic responses) associated with metal debris from implanted metallic devices may be due in part to metal ions released from implants or from particulate debris. Copyright 2002 Wiley Periodicals, Inc.

PMID: 11920666 [PubMed - indexed for MEDLINE]







	ACRI	ı	up	440		of Med	licine 🖎	L.IVI	
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	ОМІМ	PMC	alamudi,	Books
Search P	ubMed	for	Preview	······································	History	Cli	<b>Go €Cl</b> ipboard	ear	Details
Aboul Entrez	*5	Limits		Alidex	<u> </u>	2018			arr of
		<b>Display</b> Abstra	∍ct	Show	w: 20 Sort	<u> </u>	Send to	Text	
Text Version									

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NI M Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ 1: J Bone Joint Surg Br. 2004 Jan;86(1):27-30.

Related Articles, Links

### Pain in the well-fixed, aseptic titanium hip replacement. The role of corrosion.

Hallam P, Haddad F, Cobb J.

University College London Hospitals, England, UK.

We have investigated nine patients with cemented Furlong (JRI, London, UK) titanium hip replacements who presented with early pain despite a well-fixed, aseptic prosthesis. All were followed up clinically and radiologically at regular intervals. Pain was located in the thigh and was worse at night. Radiographs showed cortical hypertrophy of the femur around the tip of the stem. Eight of the nine patients subsequently required single-stage revision using an uncemented prosthesis, which relieved the pain. At revision, the pH of the tip of the stem was found to be highly acidic with macroscopic evidence of corrosion consisting of multiple layers of titanium oxides when studied by Xray dispersive analysis. Cemented titanium implants have a potential for crevice corrosion leading to cortical hypertrophy and intractable pain.

PMID: 14765860 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort Send to Text

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer







	4	an The		or Medicine	15111
Entrez PubMed	Nucleotide	Protein Genome	Structure	OMIM PMC	Journals Books
Search PubMed	for			Go C	ear
Searcit	Limits	Preview/Index	History	Clipboard	Details
About Entrez	Display Abstra	ct	Show: 20 Sor	t 🚽 Send to	Text 💽
Tout Vorsion					

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents **NLM Gateway** TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ **1:** J Arthroplasty. 2002 Oct;17(7):893-5.

Related Articles, Links

ELSEVIER FINE TEXT AND INC.

Serum levels of cobalt and chromium in a complex modular total hip arthroplasty system.

Harding I, Bonomo A, Crawford R, Psychoyios V, Delves T, Murray D, McLardy-Smith P.

Nuffield Orthopaedic Centre, Headington, Oxford, United Kingdom. ianjharding@hotmail.com

There is concern that modularity in a total hip arthroplasty system increases serum cobalt and chromium ion levels. This study measures the serum cobalt and chromium levels in patients with an Oxford Universal Hip (Corin, Cirencester, United Kingdom), which has a modular sliding mechanism; patients with a similarly manufactured hip with no sliding mechanism; and a control group. Loosening was excluded clinically and radiologically. Arthroplasty patients had statistically higher levels of serum cobalt and chromium than controls, but there was no significant difference in levels between the implanted groups. Copyright 2002, Elsevier Science (USA)

PMID: 12375249 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort	Send to Text
Control of	SHOW. I	

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer

81	<b>VCBI</b>	ļ	Pub	Me	d	Na L of Me	itional ibrary dicine	TIM	
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	OMIM	PMC	Journals	Books
Search Pu	ıbMed	for					Go C		
Jearen J		Limits	Prev	iew/Index	History	С	lipboard	De	tails
About Entrez	<u></u>	Display Abst	ract	Show	w: 20 📓 So	rt Ş	Send to	Text	
Tout Varsian									

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ 1: J Trace Elem Electrolytes Health Dis. 1992 Dec;6(4):239-43.

Related Articles, Links

# Nickel-, chrom- and cobalt-concentrations in human tissue and body fluids of hip prosthesis patients.

Hennig FF, Raithel HJ, Schaller KH, Dohler JR.

Department of Accidental Surgery, University of Erlangen-Nurnberg, Fed. Rep. of Germany.

The world-wide experience with millions of metallic implantations suggests the biocompatibility of modern alloys, commonly made of cobalt, chromium and nickel. There is, however, little information available on the internal metal exposure resulting from implants. In this study we assessed the metal concentrations in body fluids and tissue samples (muscle, bone) of patients who had undergone total hip replacement. Our patients were divided up into two groups. One group had firmly fixed implants two years after surgery. The other group had loose implants of the same Co-Cr-Mo alloy. Urine analyses revealed an increased renal elimination of nickel, chromium and cobalt. Cobalt and nickel exceeded the upper normal value. In serum the concentrations of nickel and chromium were normal or slightly elevated, the cobalt concentrations were significantly elevated. In some cases tissues adjacent to the implant showed extremely high concentrations of chromium and cobalt. This finding was also obtained in tissues that had no direct contact with the arthroplasty. The findings suggest that alloys of prostheses can undergo corrosion and release metal ions.

PMID: 1304233 [PubMed - indexed for MEDLINE]


Display Abstract	Show: 20 Sort Send to Text

Write to the Help Desk NCBI | NLM | NiH Department of Health & Human Services Freedom of Information Act | Disclaimer

Apr 19 2004 06:53:45

TILLED CLEAN TONNESS

Privacy Policy



Publication Types:

Case Reports

PMID: 9061164 [PubMed - indexed for MEDLINE]

operation (P < 0.03). The mean survival time was 13.8 years.

Display Abstract Show: 20 Sort Send to Text

operation and was significantly less (11.9 years) in those 61 years old or older at

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer

Books







PMC

Journals

OMIM

Entrez	ubMed	for			Go Clear	
Search L	-	Limits	Preview/Index	History	Clipboard	Details
About Entrez		Display Abstract		Show: 20 🖫 Sort	Send to To	ext 🔄

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ 1: Clin Orthop. 1996 Aug;(329 Suppl):S99-105.

Related Articles, Links

## Contemporary total hip replacement with metal on metal articulation.

Structure

Hilton KR, Dorr LD, Wan Z, McPherson EJ.

Center for Arthritis and Joint Implant Surgery, USC University Hospital, Los Angeles, CA, USA.

Between 1991 and 1994, 74 patients received total hip replacements with metal on metal articulation. The results of these patients with 74 hips who had a 6-month to 4-year (average, 2.2 years average) followup are reported. Patients were prospectively evaluated by the Harris hip score, a patient self assessment form, and radiographs. The average postoperative Harris hip score was 91. Patient self assessment forms showed that 95% of the patients scored their results as excellent or good. No patient had revision for loosening, but 1 underwent revision surgery for recurrent dislocation. Serial radiographs have not revealed loosening or osteolysis. Wear could not be measured radiographically. Twenty-seven of the patients had bilateral total hip replacements with 1 hip being metal on polyethylene; the patients could not detect any difference between the 2 hips. The satisfactory short term results from the contemporary metal on metal articulation investigated in this study are encouraging and warrant continued study.

PMID: 8769327 [PubMed - indexed for MEDLINE]

Display Abstract Show: 20 Sort Send to Text	Display Abstract	
---	------------------	--

Write to the Help Desk

NCBI | NLM | NiH

Department of Health & Human Services

Freedom of Information Act | Disclaimer







Entrez	PubMed '	Nucleotide	Protein	Genome	Structure	MIMO	PMC	Journals	Bı
Search	PubMed	for Limits	Previ	ew/index	History	CI	Go C pboard	l <b>ear</b> De	etails
About Entre	ez	Display Abstr	act	Show	w: 20 🔻 Sor	CONTROL AND	Send to	Text	
Text Version	n	□ 1: Clin Orth	op. 1996 .	Aug;(329 Sı	ippl):S217-32	2.	Rela	ated Articles	, Links

Entrez PubMed Overview Heip | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services
Journals Database
MeSH Database
Single Citation Matcher
Batch Citation Matcher
Clinical Queries
LinkOut
Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

Biologic effects of cobalt chrome in cell and animal models.

Howie DW, Rogers SD, McGee MA, Haynes DR.

Department of Orthopaedics and Trauma, Royal Adelaide Hospital, Australia.

The literature on animal and cellular models used to study the response to cobalt chrome alloy implants and wear and corrosion products is reviewed. Animal studies show that in solid form cobalt chrome alloy is relatively well tolerated. Injections of large numbers of particles in a single bolus lead to acute inflammation and necrosis, followed by a chronic inflammatory response. Macrophages are the predominant cell type and may persist in the tissues for years. Long term studies have failed to confirm the induction of tumors. In vitro studies confirm the toxic effects of cobalt chrome alloy corrosion products and wear particles, especially cobalt, and show that intracellular corrosion is an important mechanism for early release of cobalt ions. In vitro studies show that cobalt chrome alloy particles induce the release of inflammatory mediators from macrophages before causing cell death. These mediators have significant effects on osteoblastlike cells, as well as inducing bone resorption. Variations in the cell types, implantation site, and characteristics of the particles used in experimental models make interpretation of the results difficult. Standardized methods to control for size, shape, and number of particles for testing are proposed. It is important that in vitro and in vivo findings not be taken in isolation, but be compared with the results of human studies.

PMID: 8769336 [PubMed - indexed for MEDLINE]

			*
OBS Y WIND SPEED GA	[00 ] [00 d	Sond to Te	· -
Display Abstract	Show: 20 Sort	Send to Te	TX:

Write to the Help Desk
NCB| | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer







	***			•		
Entrez	PubMed	Nucleotide	Protein Genom	e Structure	OMIM PMC	Journals Books
Search F	ubMed	for			Go	Clear
Jearen	MANAGERY 2	Limits	Preview/Index	History	Clipboard	Details
About Entrez		<b>Display</b> Abstra	ct	Show: 20 Si	ort Send	Text

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Afeits ClinicalTrials.gov PubMed Central

Privacy Policy

☐ 1: J South Orthop Assoc. 2003 Summer;12(2):106-11.

Related Articles, Links

## Alternative bearing surfaces in total hip arthroplasty.

Inzerillo VC, Garino JP.

University of Pennsylvania School of Medicine, Philadelphia, PA, USA.

Polyethylene wear and extension of indications of total hip arthroplasty into younger and younger age groups have pushed manufacturers to develop more durable bearing surfaces. Standard polyethylene, the plastic used for the first 3 decades of hip replacement, virtually ceases to exist in its original form. Modifications of the processing, including sterlization in an inert environment and cross-linking, have demonstrated some improvements in wear. Hard-on-hard bearings such as ceramic-onceramic and metal-on-metal also have demonstrated extremely low wear. This article reviews the pros and cons of the alternative bearing options available to assist in the proper bearing selection for a particular patient.

PMID: 12882250 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort Send to Text

Write to the Help Desk

NCBI | NLM | NIH

Department of Health & Human Services
Freedom of Information Act | Disclaimer

Apr 19 2004 06:53:45

304







Sant Comment		_				OZ 1-LC			
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	ОМІМ	PMC	alamuol,	Воокѕ
Search P	ubMed	for _	an orași	······································			Go C	lear	
50010.1 m		Limits	Previe	w/Index	History	CI	lipboard		Details
About Entrez		<b>Display</b> Abstr	act	Shov	<sub>w:</sub> 20 🖫 Sor	i Ž	Send to	Text	
Text Version									

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

**PubMed Services** Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ 1: J Bone Joint Surg Am. 1991 Dec;73(10):1475-86.

Related Articles, Links

Comment in:

J Bone Joint Surg Am. 1992 Oct;74(9):1431-2.

Release and excretion of metal in patients who have a total hipreplacement component made of titanium-base alloy.

Jacobs JJ, Skipor AK, Black J, Urban R, Galante JO.

Department of Orthopedic Surgery, Rush-Presbyterian-St. Luke's Medical Center, Chicago, Illinois 60612.

Serum concentration and urinary excretion of titanium, aluminum, and vanadium were measured for patients who had a well functioning cementless primary total hip replacement of one of two different designs, for patients who had a loose total hip replacement that was to be revised, and for control subjects who had no implant. Serum concentrations of titanium were elevated approximately twofold in the patients who had a loose implant, compared with the values for the control subjects. No major differences in terms of urine concentration of titanium, serum concentration of aluminum, or urine concentration of aluminum were observed among any of the groups that were studied. Concentrations of vanadium were uniformly low in all groups.

PMID: 1844129 [PubMed - indexed for MEDLINE]

And the state of t	
St. American Care March 1990 Company	A Control of Property of Control
Display Abstract	Show: 20 Sort Send to Text

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer







				•					
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	OMIM	PMC	Journals,	Book
Search P	ubMed	for			A CONTRACTOR OF THE PROPERTY O		Go C	lear.	
		Limits	Previe	w/Index	History	С	lipboard	Deta	ails
About Entrez		<b>Display</b> Abstrac	ct	Sho	w: 20 🔀 Soi	t X	Serid to	Text	M

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

**PubMed Services** Journals Database MeSH Database Single Citation Matcher **Batch Citation Matcher** Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ 1: Clin Orthop. 1996 Aug;(329 Suppl):S256-63.

Related Articles, Links

Cobalt and chromium concentrations in patients with metal on metal total hip replacements.

Jacobs JJ, Skipor AK, Doorn PF, Campbell P, Schmalzried TP, Black J, Amstutz HC.

Department of Orthopedic Surgery, Rush Arthritis and Orthopedics Institute, Chicago, IL, USA.

There has been a resurgence of interest in the use of metal on metal bearings in total hip arthroplasty. Although the use of metal on metal bearing couples would eliminate or substantially reduce particulate polyethylene generation (depending on the presence or absence of polyethylene in the implant system), there is concern about the potential for increased particulate and ionic metal generation in comparison with polyethylene on metal bearings. These metallic degradation products may be transported away from the implant site and distributed systemically. Chromium concentrations in the serum and urine and cobalt concentrations in the serum were measured in subjects with cobalt chromium alloy metal on metal total hip replacements and in controls without implants. Eight subjects with long term (> 20 years) McKee-Farrar total hip replacements had 9fold elevations in serum chromium, 35-fold elevations in urine chromium, and at least 3fold elevations in serum cobalt concentrations in comparison with controls. Six subjects with short term (< 2 years) metal on metal surface replacement arthroplasties had 3fold elevations in serum chromium, 4-fold elevations in urine chromium, and 4-fold elevations in serum cobalt concentrations in comparison with subjects with McKee-Farrar implants. Although the toxicologic importance of these trace metal elevations has not been established, serum and urine metal concentrations may be useful markers for the tribologic performance of metal on metal bearings.

PMID: 8769339 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort Send to Text

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer







Entrez	PubMed	Nucleotide	Protein	Genome	Structure	MIMO	PMC	Journals	RDOW
<del></del>	PubMed	for					Go C	lear	
Search	J	Limits	Previev	v/Index	History	Ci	ipboard	Det	tails
About Entr	ez	<b>Display</b> Abstrac	ct	Sho	w: 20 So	rt 🎇	Send t	o⊭ Text	¥

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ 1: J Bone Joint Surg Am. 1998 Oct;80(10):1447-58.

Related Articles, Links

Full text article at www.ejbjs.org

Metal release in patients who have had a primary total hip arthroplasty. A prospective, controlled, longitudinal study.

Jacobs JJ, Skipor AK, Patterson LM, Hallab NJ, Paprosky WG, Black J, Galante

Department of Orthopedic Surgery, Rush Arthritis and Orthopedics Institute, Rush Medical College, Rush-Presbyterian-St. Luke's Medical Center, Chicago, Illinois 60612, USA. jacobs@ortho4.pro.rpslmc.edu

There is an increasing recognition that, in the long term, total joint replacement may be associated with adverse local and remote tissue responses that are mediated by the degradation products of prosthetic materials. Particular interest has centered on the metal-degradation products of total joint replacements because of the known toxicities of the metal elements that make up the alloys used in the implants. We measured the concentrations of titanium, aluminum, cobalt, and chromium in the serum and the concentration of chromium in the urine of seventy-five patients during a three-year prospective, longitudinal study. Twenty patients had had a so-called hybrid total hip replacement (insertion of a modular cobalt-alloy femoral stem and head with cement and a titanium acetabular cup without cement), fifteen had had insertion of an extensively porous-coated cobalt-alloy stem with a cobalt-alloy head and a titaniumalloy socket without cement, and twenty had had insertion of a proximally porouscoated titanium-alloy stem with a cobalt-alloy head and a titanium socket without cement. The remaining twenty patients did not have an implant and served as controls. The results of our study showed that, thirty-six months postoperatively, patients who have a well functioning prosthesis with components containing titanium have as much as a threefold increase in the concentration of titanium in the serum and those who have a well functioning prosthesis with cobalt-alloy components have as much as a fivefold and an eightfold increase in the concentrations of chromium in the serum and urine, respectively. The predominant source of the disseminated chromium-degradation products is probably the modular head-neck junction and may be a function of the geometry of the coupling. Passive dissolution of extensively porous-coated cobalt-alloy stems was not found to be a dominant mode of metal release. CLINICAL RELEVANCE: Increased concentrations of circulating metal-degradation products derived from orthopaedic implants may have deleterious biological effects over the long term that warrant investigation. This is a particularly timely concern because of recent clinical trends, including the reintroduction of metal-on-metal bearing surfaces and the increasing popularity of extensively porous-coated devices with large surface areas of exposed metal. Accurate monitoring of the concentrations of metal in the serum and urine after total hip replacement also can provide insights into the mechanisms of metal release. Our findings suggest that fretting corrosion at the head-neck coupling is an important source of metal release that can lead to increased concentrations of chromium in the serum. Determinations of the concentrations of metal in the serum and urine may be useful in the diagnosis of patients who are symptomatic after a total joint replacement as increased levels are indicative of at least one mode of mechanical dysfunction (for example, fretting corrosion) of the device.


PMID: 9801213 [PubMed - indexed for MEDLINE]

= <del></del>		
Display Abstract	Show: 20 Sort	Send to Text

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer

PubMed Central

Privacy Policy



were present in 52% of the surviving prostheses. Clinical scores showed weak correlation with the radiographic loosening in both groups, and 18 McKee-Farrar and 8 Charnley prostheses were still considered satisfactory by the patients. The mean annual linear polyethylene wear was 0.12 mm. Osteolytic lesions were observed in association with 2 McKee-Farrar and 5 Charnley prostheses in surviving hips. The long term results of the McKee-Farrar prosthesis are comparable with those of the low friction arthroplasty in this series. Wear of the polyethylene bearing and accumulation of polyethylene particles in the periprosthetic tissue may become an increasing problem. Second generation all metal implants seem to be worth considering in patients with long

PMID: 8769323 [PubMed - indexed for MEDLINE]

life expectancy.

and the company of the contract of the contrac	
	Beauty and the second second
Display Abstract	Show: 20 Sort Send to Text

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer



Entrez PubMed

New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database

Single Citation Matcher

Batch Citation Matcher

Related Resources

Order Documents

Consumer Health

ClinicalTrials.gov PubMed Central

NLM Gateway TOXNET

Clinical Alerts

Privacy Policy

Clinical Queries LinkOut

Cubby

Overview Help | FAQ Tutorial





	MORL		u	141 C	<u>u</u>	of Me	dicine E	LAYE	
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	OMIM	PMC	Journals	В
Search F	PubMed	for Limits	Previ	ew/Index	History	Cli	∫ <b>∘Go</b> ∐⊪ <b>C</b> pboard	l <b>ear</b> Det	tails
About Entre	∋z	<b>Display</b> Abstr	act	Show	w: 20	t Santanian na samaha santa	Send to	Text	ужинсо.
Text Versio	n	□1: Am J Orth	iop. 1998	3 Apr;27(4):2	283-92.		Rela	ated Articles,	Links

Hard bearing surfaces in total hip arthroplasty.

Jazrawi LM, Kummer FJ, Di Cesare PE.

Department of Orthopaedic Surgery, Hospital for Joint Diseases Orthopaedic Institute, New York, New York, USA.

Periprosthetic osteolysis and aseptic loosening are serious problems affecting the outcome of total joint replacement. Polyethylene particulate debris generated from metal-on-polyethylene bearing surfaces and the resulting biologic response to this debris are thought to be largely responsible. As a result, there has been a renewal of interest in hard bearing surfaces for total joint arthroplasty, including both metal-on-metal and ceramic-on-ceramic components. The new-generation all-ceramic and all-metal prostheses have demonstrated, both clinically and in the laboratory, lower friction and wear rates than metal-on-polyethylene bearing surfaces. Theoretically, lower wear rates result in less particulate debris and decreased inflammatory response. Despite excellent tribologic (lubrication, friction, wear) properties, metal-on-metal bearings raise associated issues of metal sensitivity and toxicity. For ceramic-on-ceramic bearing surfaces, issues of ceramic quality and the possibility of brittle fracture must be considered.


Publication Types:

- Review
- Review, Tutorial

PMID: 9586727 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort	Send to Text

Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer



Related Resources
Order Documents
NLM Gateway
TOXNET
Consumer Health
Clinical Alerts
ClinicalTrials.gov
PubMed Central

Batch Citation Matcher Clinical Queries

LinkOut

Cubby

Privacy Policy

an interest in developing polyethylene with improved wear characteristics, as well as a renewed interest in alternative bearing surfaces for total joint arthroplasty, including ceramic-polyethylene, metal-metal, and ceramic-ceramic articulations. These alternative surfaces have demonstrated less friction and lower wear rates than metal-on-polyethylene bearing surfaces in both clinical and laboratory experiments. Clinical results, although only short- to mid-term, have been encouraging. Alternative bearing surfaces, with lower wear rates and less particulate debris formation, may have the potential to improve total joint arthroplasty survivorship by decreasing periprosthetic osteolysis, especially in younger, high-demand patients.

The biologic response to polyethylene particulate debris generated from metal-on-

polyethylene bearing surfaces is thought to be largely responsible for periprosthetic

osteolysis and aseptic loosening in total joint arthroplasty. As a result, there has been

#### Publication Types:

- Review
- · Review, Tutorial

PMID: 9682082 [PubMed - indexed for MEDLINE]

Application for the control of	201 201	mediane from the Asian
Display Abstract	Show: 20 Sort	Send to Text


Write to the Help Desk
NCBI | NLM | NiH
Department of Health & Human Services
Freedom of Information Act | Disclaimer

S NCBI	P	ubline	ed	National Library of Medicine	
Entrez PubMed	Nucleotide	Protein Genome	Structure	OMIM PMC	Journals Books
DubMed	for			GoClea	r j
Search I Labined	Limits	Preview/Index	History	Clipboard	Details
About Entrez	Display Abstrac	st Sh	ow: 20 Sort	✓ Send to	Text
Text Version					
Entrez PubMed Overview Help   FAQ Tutorial New/Noteworthy E-Utilities	Springer La	nk Mikika eller with meta	al-on-metal a	Epub 2002 Dec 19.  rticulation: clin s of short-term	Related Articles, Links ical, results.
PubMed Services Journals Database MeSH Database		P, Petsinis G, Re			
Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut	korovess@o	tenet.gr		Andreas' Patras, Gre	
Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central  Privacy Policy	patients, where metal-on-metal-on-metal-on-metal-on-mevaluated confusion for their lass. The mean factors was apatients we aseptic loose Dislocation Periarticula IV), but wit could be idevidence for Zweymus 96.8%, resuggested to set of the se	no received 350 thirdetal articulation for linically, radiological was 55+/-9 years, (in the follow-up evaluation of the prosthesis (Mirra gueller screw socket pectively. CONCLUST of the prosthesis of a converse the prosthesis of the prosthesis o	d-generation Zwe primary or second ly and histologica range 25-70 years on for reasons unconths (range 37-9) ed to 96+/-4 postoved postoperative satisfied with the series. Revision wor technical error ecurred in the earloserved in 30% of ability. During revy formed hip joint grades 1 and 2) wand stem 7.6 years 10NS: This study asty with metal-or entional polyethyles.	HOD:. A total of 266 cymueller-SL total hip dary osteoarthritis, wally. The age of the particular particular of the particular of the particular of the operatively. The invalue of the operatively (p<0.001). In all result of the operative particular of the operative period of the hips (5% Brook ision surgery, no man membrane; howeve as seen in all revised as after implantation showed that the showed on ceramic articular operations of the on ceramic articular operations of the operation of the on ceramic articular operations of the operati	as followed and atients at the time ents did not return peration. RESULTS perative Harris hip alidity of the peration. There was no hips because of implantation. In additional and croscopic metalosis of microscopic metalosis of microscopic lips. The survival was 99.4% and port-term results of were just as ulation, while the
	complication PMID: 125	82789 [PubMed - in	dexed for MEDLIN	NE]	

Send to Text

Display Abstract

Show: 20 Sort



PMID: 1847099 [PubMed - indexed for MEDLINE]


synovium.

Display Abstract	Sh	ow: 20	Sort	3	00114 10	Text	-]


macrophage inflammatory activity is down-regulated in periprosthetic

production in response to stimulation between primary and revision groups was statistically significant. The observation of a chronic moderate level of activation and the lack of responsiveness to a potent stimulator suggests that

> Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer







			7				
i	Nucleotide	Protein Genoi	ne Structure	MIMO	PMC	Journals B	ook
	for		MINOR INC.		Go Cle	el.	
	Limits	Preview/Index	: History	С	lipboard	Details	
<u>)</u>	<b>Display</b> Abstrac	t 🚼	Show: 20	Sort 🗲	Send to	Text	

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

**PubMed Services** Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ 1: J Bone Joint Surg Br. 1996 Jan;78(1):18-21.

Related Articles, Links

### Systemic release of cobalt and chromium after uncemented total hip replacement.

Kreibich DN, Moran CG, Delves HT, Owen TD, Pinder IM.


Freeman Hospital, Newcastle upon Tyne and Southampton General Hospital, England.

We measured the levels of cobalt and chromium in the serum in three groups of patients after uncemented porous-coated arthroplasty. Group 1 consisted of 14 consecutive patients undergoing revision for aseptic loosening. Group 2 comprised 14 matched patients in whom the arthroplasty was stable and group 3 was 14 similarly matched patients with arthritis awaiting hip replacement. Specimens were analysed using atomic absorption spectrophotometry. Aseptic loosening of a component resulted in a significant elevation of serum cobalt (p < 0.05), but not of serum chromium. The relative risk of a component being loose, if the patient had a serum cobalt greater than 9.0 nmol/l, was 2.8.

PMID: 8898120 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort	Send to Text
200 (Contracting Contracting C	SHOW.	

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer



Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer

Sort

Show: 20


Abstract

Display

Jun 1 2004 06:47:16

Text

Send to



☐ **1:** J Orthop Res. 1997 Jan;15(1):40-9.

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy F-Utilities

**PubMed Services** Journals Database MeSH Database Single Citation Matcher **Batch Citation Matcher** Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

Human monocyte/macrophage response to cobalt-chromium corrosion products and titanium particles in patients with total joint replacements.


Lee SH, Brennan FR, Jacobs JJ, Urban RM, Ragasa DR, Glant TT.

Department of Orthopedic Surgery, Rush Arthritis and Orthopedics Institute, Rush Medical College, Rush-Presbyterian-St. Luke's Medical Center, Chicago, IL 60612, USA.

The responses of human peripheral blood monocytes of 10 normal volunteers and 14 patients with total hip replacements to particles of commercially pure titanium and chromium orthophosphate (a corrosion product from cobalt-chromium alloy implants) were studied. In addition, these phagocytosable particles were added to cultured mononuclear cells isolated from the interfacial membrane of 14 patients with failed implants. Peripheral blood monocytes from patients who had had a total hip replacement produced significantly higher levels of interleukin-1 (both interleukin-1 alpha and interleukin-1 beta) and prostaglandin E2 following particulate stimulation than those from normal volunteers. Supernatants from both titanium and chromium orthophosphate-stimulated peripheral blood monocytes from the volunteers and patients with total hip replacement induced bone resorption (assayed in organ cultures of newborn mouse calvariae) and the proliferation of human fibroblasts. The levels of bone resorption were significantly higher (p < 0.05) in patients with implants than in normal volunteers. There were no significant differences in the responses of cells between patients with focal osteolysis and those without osteolysis. Interfacial membrane mononuclear cells also produced high levels of interleukin-1 alpha, interleukin-1 beta, and prostaglandin E2 and expressed bone resorptive activities following stimulation with either titanium or chromium orthophosphate. More importantly, interfacial membrane mononuclear cells "spontaneously" produced high levels of prostaglandin E2 that were comparable with the response of peripheral blood monocytes stimulated with particulate wear debris. The clinical relevance of this study is 2-fold. First, mononuclear cells from patients with total hip replacement were some-how "sensitized" to metal particles in comparison with mononuclear cells from individuals without an implant. Second, the chromium orthophosphate corrosion product was a potent macrophage/monocyte activator and may contribute to macrophage-mediated osteolysis and aseptic loosening.

PMID: 9066525 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort Send to Text	
	I Snow: I	_



Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer

20

Abstract

Display

Sort

Apr 19 2004 06:53:45

Send to

Text

Entrez PubMed Overview Help | FAQ

New/Noteworthy E-Utilities

PubMed Services

Journals Database MeSH Database Single Citation Matcher

Batch Citation Matcher Clinical Queries

Related Resources Order Documents

NLM Gateway

Clinical Alerts

ClinicalTrials.gov

PubMed Central

Privacy Policy

TOXNET Consumer Health


IsinotuT

1 inkOut

Cubby







				•	•				
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	MIMO	PMC	Journals	Books
Search	PubMed	for	· · · · · · · · · · · · · · · · · · ·				Go C	lear	
		Limits	Prev	/iew/Index	History	Cli	Clipboard		Details
About Entre	ez	<b>Display</b> Abs	stract		Show: 20 So	rt 🖫	Send to	Text	Š
Text Version	n								

☐ **1:** Arch Orthop Trauma Surg. 1986;105(5):263-7.

Related Articles, Links

Development of heterotopic ossification around the hip. A longterm follow-up of patients who underwent surgery with two different types of endoprostheses.


Lindholm TS, Viljakka T, Vankka E, Popov L, Lindholm TC.

Heterotopic ossification has been reported in many pathological situations, most important clinically as a sequel to hip arthroplasty and spinal trauma. The etiology of heterotopic ossification is yet not clear, but the disease is supposed to be connected with trauma. Heterotopic bone was found in 53% (1.2% with the severe form) of 623 patients operated on at the Orthopaedic Hospital of the Invalid Foundation, Helsinki, Finland; the operations included 849 arthroplasties. The rate of heterotopic ossification was higher after revision arthroplasty, following operation of the contralateral side, in men, and in primary coxarthrosis, and the incidence was higher with the Brunswik (metal-on-plastic) endoprosthesis than in the McKee-Farrar type (metal-on-metal). Heterotopic bone formation generally seemed to increase and to be more manifest during long-term observation.

PMID: 3778160 [PubMed - indexed for MEDLINE]

Display: Abstract	Show: 20	Sort 🔻	Send to Text
34 (A.A.) (1.10 (A.A.)	311074	B	

Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer




PMID: 11742463 [PubMed - indexed for MEDLINE]

Display Abstract Show: 20 Sort Send to Text

Write to the Help Desk
NCBI | NLM | NiH
Department of Health & Human Services
Freedom of Information Act | Disclaimer

Privacy Policy



Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer

Jun T 2004 06:47:16

320





National Library of Medicine

	A COURT	1		MA-	<u>~</u>	of Me	dicine 🕰	Live	
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	OMIMO	PMC	Journals	Books
Search Pu	ıbMed	for					Go C	lear	
JCG1 CII Remoun		Limits	Previe	w/Index	History	C	lipboard	Deta	ils
About Entrez		Display Abstra	act	Shov	v: 20 × So	n 📝	Send to	Text	

Text Version

Entrez PubMed Overview Help [ FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

**1:** Clin Orthop. 2003 Jan; (406):282-96.

Related Articles, Links



Metal-on-metal versus polyethylene in hip arthroplasty: a randomized clinical trial.

MacDonald SJ, McCalden RW, Chess DG, Bourne RB, Rorabeck CH, Cleland D, Leung F.

Division of Orthopaedic Surgery, University of Western Ontario & London Health Sciences Centre (University Campus) 339 Windermere Road, London, Ontario, N6A 5A5, Canada. steven.macdonald@lhsc.on.ca


A prospective, randomized, blinded clinical trial was done to evaluate polyethylene versus metal bearing surfaces in total hip replacement. Forty-one patients were randomized to receive either a metal (23 patients) or a polyethylene (18 patients) insert. The femoral and acetabular components were identical with the acetabular insert the only variable. Patients were assessed preoperatively and postoperatively using radiographs, multiple outcome measures (Western Ontario MacMaster University Score, Harris hip score, Short Form-12), erythrocyte metal ion analysis (cobalt, chromium, titanium), and urine metal ion analysis (cobalt, chromium, titanium). Patients were followed up for a minimum of 2 years (mean 3.2 years; range, 2.2-3.9 years). There were no differences in radiographic outcomes or outcome measurement tools between patients. Patients receiving a metal-on-metal articulation had significantly elevated erythrocyte and urine metal ions compared with patients receiving a polyethylene insert. Patients who had metal-on-metal inserts had on average a 7.9-fold increase in erythrocyte cobalt, a 2.3-fold increase in erythrocyte chromium, a 1.7-fold increase in erythrocyte titanium, a 35.1-fold increase in urine cobalt, a 17.4-fold increase in urine chromium, and a 2.6-fold increase in urine titanium at 2 years followup. Patients receiving a polyethylene insert had no change in erythrocyte titanium, urine cobalt, or urine chromium and a 1.5-fold increase in erythrocyte cobalt, a 2.2-fold increase in erythrocyte chromium, and a 4.2-fold increase in urine titanium. Forty-one percent of patients receiving metal-on-metal articulations had increasing metal ion levels at the latest followup.

#### Publication Types:

- Clinical Trial
- Randomized Controlled Trial

PMID: 12579029 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort	Send to Text →
------------------	---------------	----------------



Write to the Help Desk

NCBI | NLM | NIH

Department of Health & Human Services

Freedom of Information Act | Disclaimer

Sort

Show: 20

Display

Abstract

Jun 1 2004 06:47:16


Text

Send to

322







nal try ne NLM	

Entrez	PubMed	Nucleotide	Protein Genom	e Structure	OMIM PIV	C Journals	Boo
Search	PubMed	for	www.s.c.iii.		<u>G</u> e	Clear	
DCG(C)		Limits	Preview/Index	History	Clipboa	rd	Details
bout Entr	ez	Display Abstrac	et 🔄	Show: 20 So	ort 💆 Se	nd to Text	Š

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher **Batch Citation Matcher** Clinical Queries LinkOut Cubby

Related Resources Order Documents **NLM Gateway** TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

1: Clin Orthop. 1996 Aug;(329 Suppl):S128-40.

Related Articles, Links

In vivo wear of three types of metal on metal hip prostheses during two decades of use.

McKellop H, Park SH, Chiesa R, Doorn P, Lu B, Normand P, Grigoris P, Amstutz

J. Vernon Luck Orthopaedic Research Center at the Orthopaedic Hospital, Los Angeles, CA, USA.

Wear was analyzed on 21 metal on metal hip replacements, including McKee-Farrar, Muller, and Ring, that were retrieved from patients after as many as 25 years. Light and scanning electron microscopy indicated that early wear included substantial third body abrasion, possibly from particles generated while scratches from the original polishing were being eradicated and from dislodged surface carbides. However, the main contact zones were eventually worn smoother than the original surfaces. Wear was quantified by digitizing the shapes of the components on a coordinate measuring machine and identifying those areas that deviated from the original spheric surface. On the femoral heads, wear was typically concentrated in the superomedial region, that is, on the load axis. Three cases also had substantial wear inferiorly, but there were no cases with circumferential (equatorial) wear. The long term wear rates averaged approximately 6 micrometers per year or less and produced an average of approximately 6 mm3 of metallic wear debris per year or less. Wear rate tended to increase as clearance increased over the range of 127 to 386 micrometers, and a McKee-Farrar prosthesis with the extreme clearance of 1.7 mm wore approximately 16 times faster than the average, but there was no apparent relationship between clearance and time to revision. Larger McKee-Farrar balls had less volumetric wear, on average, than smaller balls, and the Muller balls had the greatest wear, which may have been due to contact with the edges of recesses machined into the bearing zones of the Muller cups.

PMID: 8769330 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort	Send to Text

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer

Text Version

Overview

Help | FAQ Tutorial

New/Noteworthy E-Utilities

PubMed Services Journals Database

Batch Citation Matcher Clinical Queries

Related Resources

Order Documents

Consumer Health

PubMed Central

Privacy Policy

NLM Gateway TOXNET

Clinical Alerts ClinicalTrials.gov

MeSH Database Single Citation Matcher

LinkOut

Cubby







				•					
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	OMIM	PMC	Journals	Вро
Search P	ubMed	for			voca na voca n		Go 🖰	Clear	
Jedien III.		Limits	Previev	v/Index	History	Clip	board	[	Details
About Entrez		<b>Display</b> Abstra	ct	Sho	w: 20 👺 Sori		Send	to Text	is.

Related Articles, Links

## Metal on metal surface replacement of the hip. Experience of the McMinn prothesis.

McMinn D, Treacy R, Lin K, Pynsent P.

Midland International Orthopaedic Service, Birmingham Nuffield Hospital, United Kingdom.

The historical failure of surface replacement has been due to the production of wear debris with subsequent bone resorption, loosening, and failure. To avoid these problems, a surface replacement using a metal on metal bearing allowing thin components and femoral design and instrumentation to avoid varus alignment has been designed. Two hundred thirty-five joints have been resurfaced with this prosthesis in almost 5 years. There have been no femoral neck fractures and no dislocations. There have been 4 designs differing in the method of fixation. In the press fit group, 6 of 70 hips had to be revised for aseptic loosening. In the cemented group, debonding of the cup occurred in 3 of 43 cases. Six patients had hydroxyapatite coated components and have had excellent clinical outcomes. The current design uses a peripherally expanded hydroxyapatite coated cup and a cemented metal head; 116 of this design have been implanted during a 19-month period with excellent outcome. Despite short followup the authors are hopeful that the combination of a polar metal on metal bearing with appropriate fixation will yield a method of preserving bone stock in the younger patient requiring arthroplasty.

PMID: 8769326 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort	Send to Text

Write to the Help Desk

NCBI | NLM | NIH

Department of Health & Human Services
Freedom of Information Act | Disclaimer

Jun 1 2004 06:47:16


324

	•				
S NCBI	Pı	Mode	ed	National Library of Medicine	M
Emrez PubMed	Nucleotide Prof	ein Genome	Structure	OMIM PMC	Journals Books
Search PubMed	for		The state of the s	Go Cle	ear
JULI CII Parametria	Limits	Preview/Index	History	Clipboard	Details
About Entrez	Display Abstract	S S	how: 20 🖺 Sort	<b>y</b> Send to	Text 💆
Text Version					<u> </u>
Entrez PubMed	1: Clin Orthop. 19	96 Aug;(329 Sur	ppl):S233-43.		Related Articles, Links
Overview Help   FAQ Tutorial New/Noteworthy E-Utilities		n of cobalt o ic reactions.		ar and corrosio	on products
PubMed Services	Merritt K, Bro	own SA.			•
Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher	Radiologic Hea	alth, Office of Sci	ence and Technolo	ration, Center for I ogy, Rockville, MD,	USA.
Clinical Queries LinkOut Cubby  Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central  Privacy Policy	for the cup art years of succe articulating su is now being rand corrosion ultrahigh mole release, and be have been sho and the alloy replacements. corrosion propand corrosion replacement of at the articular polyethylene, low to be well transported from tissue and eliminations of condermatitis. Under the succession of condermatitis.	siculating with a ss. However, the rfaces continues econsidered for of the articulating cular weight policiologic fate, of nown to have limit per se should not be cobalt chroperties and continuous and he devices is review thing surfaces of the amount of not tolerated by the om the implant minated more slentoward biologic	metal head has pro- e search for improv- total hip replacements g surfaces can be yethylene couple. The considered for mium alloys are known to be evaluate ow this relates to the ed. Under the cond cobalt chromium a metallic products tre biologic system. Site and eliminated owly. The issue of ments, in ionic form	kept below that of Concern has been recorrosion and wear clating surface show recorrosion to have reaso d. The issue of cobne biologic perform dition of wear as culloys and ultrahigh ansferred to the tis Nickel and cobalt in the urine. Chrone host hypersensitivity, are known to caughypersensitivity,	lesigns for g with metallic cups e enhanced if wear the metal on raised about the r. Titanium alloys ving significant wear, otal hip onable wear and balt chromium wear ance of total hip urrently experienced molecular weight sues is sufficiently ons arc, rapidly mium is stored in the ty to these elements
	Publication Ty • Reviev • Reviev	•			
	PMID: 87693	37 [PubMed - inc	dexed for MEDLIN	<b>[</b> ]	

Send to Text

Show: 20 Sort

Display Abstract



Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer

E-Utilities

LinkOut

Cubby

PubMed Services Journals Database MeSH Database

Single Citation Matcher Batch Citation Matcher Clinical Queries

Related Resources

Order Documents

**NLM Gateway** 

Clinical Alerts

Privacy Policy

ClinicalTrials.gov PubMed Central

TOXNET Consumer Health

81	NCBI		Pub	Me	d	Li	ional brary licine	IM	
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	MIMO	PMC	Journals	Books
Search Pu	ıbMed	for					Go C	lear	
Search L.		Limits	Previ	ew/Index	History	Cli	pboard	. De	tails
About Entrez		Display Abs	tract	Sho	w: 20 🔀 Sc	nt 💆	Send to	Text	
Text Version					•				
Entrez PubM	ed	🗆 <b>1:</b> ) Natl Car	icer Inst. 19	995 Jan 4;87∙	(1):28-33.			Related A	Articles, Links
Overview Help   FAQ Tutorial New/Noteworth		Cancer	risk afte tion-bas	er hip rep	lacement v t study in S	vith met Sweden.	al imp	lants: a	. •

population-based conort study in Swed

Nyren O, McLaughlin JK, Gridley G, Ekbom A, Johnell O, Fraumeni JF Jr, Adami

Department of Cancer Epidemiology, University Hospital, Uppsala, Sweden.

BACKGROUND: Joint replacement with metal implants has been the standard procedure for surgical treatment of irreversible degeneration of hip and knee joints for more than two decades. However, reports of local malignancy after joint replacement and experimental studies that suggest a carcinogenic action of metal ions and polymethylmethacrylate (an acrylic compound used to stabilize the implant in the host) have raised concern about the possible long-term risks associated with metal implants. PURPOSE: Our aim was to study cancer risk in a Swedish cohort of patients who had hip replacement surgery during the period 1965 through 1983. METHODS: We studied the risk of cancer in a cohort of 39 154 patients (14 869 men and 24 285 women), identified in the nationwide Swedish Inpatient Register with at least one hip replacement during the period 1965 through 1983. The patients were followed through 1989 by means of record linkage to the Swedish Cancer Register. The cohort contributed a total of 327 922 person-years at risk. Standardized incidence ratios (SIRs) were computed using age-, sex-, and period-specific incidence rates derived from the entire Swedish population. RESULTS: The overall relative risk of cancer was increased by only 3%. Bone cancer--the focus of previous concerns--occurred in six cases versus 4.3 expected, and connective tissue cancer occurred in 28 cases versus 25.9 expected. Increased risks were observed for kidney cancer (SIR = 1.31; 95% confidence interval [CI] = 1.13-1.51), prostate cancer (SIR = 1.13; 95% CI = 1.04-1.22), and melanoma (SIR = 1.23; 95% CI = 1.00-1.50). The relative risk of gastric cancer steadily declined with increasing follow-up time, in both men and women ( $\overline{SIR} = 0.58$ ; 95% CI = 0.39-0.84 more than 10 years after hip replacement). CONCLUSION: In this study, the largest study to date to evaluate hip replacement and subsequent cancer risk, the overall cancer risk appears to be negligible from a public health perspective, and our results have not produced any strong evidence against the continued use of these devices. Nevertheless, the small but statistically significant increases in kidney and prostate cancers and the decrease in gastric cancer deserve further study.

PMID: 7666459 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort Send to Text	





National Library of Medicine

			•	₹				
Entrez	PubMed	Nucleotide	Protein Gener	ne Structure	MIMO	PMC	Journals	Books
Search Pub	Med	for	**************************************			Go Clea	ır.	
		Limits	Preview/Index	History	Clip	pboard	Details	i
About Entrez	<b>一</b>	Display Abstrac	et 🌠	Show: 20 S	ort 🔾	Send to	Text	MAN TO THE PROPERTY OF THE PRO

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ **1:** Radiographics. 2000 May-Jun;20(3):699-712.

Related Articles, Links

FREE full text article at radiographics.rsnajnls.org

Metal artifact reduction sequence: early clinical applications.

Olsen RV, Munk PL, Lee MJ, Janzen DL, MacKay AL, Xiang QS, Masri B.

Departments of Radiology, University of British Columbia, Vancouver General Hospital, 855 W 12th Ave, Vancouver, British Columbia, Canada.

Artifact arising from metal hardware remains a significant problem in orthopedic magnetic resonance imaging. The metal artifact reduction sequence (MARS) reduces the size and intensity of susceptibility artifacts from magnetic field distortion. The sequence, which is based on view angle tilting in combination with increased gradient strength, can be conveniently used in conjunction with any spin-echo sequence and requires no additional imaging time. In patients with persistent pain after femoral neck fracture, the MARS technique allows visualization of marrow adjacent to hip screws, thus enabling diagnosis or exclusion of avascular necrosis. Other applications in the hip include assessment of periprosthetic soft tissues after hip joint replacement surgery, postoperative assessment after resection of bone tumors and reconstruction, and localization of unopacified methyl methacrylate cement prior to hip arthroplasty revision surgery. In the knee, the MARS technique allows visualization of structures adjacent to implanted metal staples, pins, or screws. The technique can significantly improve visualization of periprosthetic bone and soft-tissue structures even in patients who have undergone total knee arthroplasty. In patients with spinal fixation hardware, the MARS technique frequently allows visualization of the vertebral bodies and spinal canal contents. The technique can be helpful after wrist fusion or screw fixation of scaphoid fractures.


PMID: 10835123 [PubMed - indexed for MEDLINE]

	20 4 5		Send to Text ¥
<b>Display</b> Abstract →	Show:  20   3	SUIL 🚰	OCHORO TEXT

Write to the Help Desk
NCBI | NLM | NiH
Department of Health & Human Services
Freedom of Information Act | Disclaimer







Entrez	PubMed	Nucleotide	Protein	Genome	Structure	OMIM	PMC	Journals 	В
Search P	ubMed	for Limits	Prev	iew/Index	History	Cli	Go C pboard	lear De	etails
About Entre	z	Display Abstra	act	Show	w: 20 Sor	t S	Send to	Text	21.73 23.43
Text Version	n	□1: J Arthropk	asty, 19!	99 Apr;14(3)	:272-80.		Rela	ated Articles	. Links

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services
Journals Database
MeSH Database
Single Citation Matcher
Batch Citation Matcher
Clinical Queries
LinkOut
Cubby

Related Resources
Order Documents
NLM Gateway
TOXNET
Consumer Health
Clinical Alerts
ClinicalTrials.gov
PubMed Central

Privacy Policy

Erratum in:

• J Arthroplasty 2000 Jan;15(1):136-7.

J Arthroplasty

Cancer incidence in Finnish hip replacement patients from 1980 to 1995: a nationwide cohort study involving 31,651 patients.

Paavolainen P, Pukkala E, Pulkkinen P, Visuri T.


National Agency for Medicines, and the Department of Surgery, The Jorvi Hospital, Espoo, Finland.

Nationwide, computer-based reporting of all arthroplasties performed in Finland was started in January 1980. Using data from these records, a cohort of 31,651 polyethylene-on-metal total hip arthroplasty (THA) patients was followed up for cancer, using Finnish Cancer Registry data, from 1980 to 1995. During follow-up, 2,367 cancers were observed. There were statistically significantly fewer cancers among the THA patients (standardized incidence ratio [SIR], 0.90; 95% confidence interval [CI], 0.87-0.93). SIRs for cancers of the lung (0.69) and stomach (0.77) were significantly below unity. There was no significantly increased risk at any site. The SIR for cancer overall in male THA patients was below unity during the first 3 years after THA but returned to unity thereafter. The low SIR among men during the first 3 years was largely because the lung cancer SIR was 0.47 (95% CI, 0.35-0.62). In women, the SIR remained around 0.93 throughout follow-up. The SIR for stomach cancer was below unity only in women (SIR, 0.67; 95% CI, 0.51-0.86). For cancer of the urinary bladder, the SIR during the first 3 years after THA was below unity but later slightly above it (SIR, 1.24 in relation to > or =3 years of follow-up; 95% CI, 0.99-1.52). For myeloma and leukemia, SIRs were greater than unity only for THA patients followed up for 3 to 9 years. The study findings, in contrast to previously reported findings, do not indicate that there is any increased risk of hematopoietic cancers after THA using polyethylene-onmetal prostheses. SIRs relating to soft tissue cancers and bone sarcomas did not differ significantly from unity. No sarcoma was observed at the site of a prosthesis. THA seems to play no major role in cancer causation.

PMID: 10220179 [PubMed - indexed for MEDLINE]

PubMed Central


Privacy Policy



Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer

ClinicalTrials.gov PubMed Central

Privacy Policy



PMID: 9322702 [PubMed - indexed for MEDLINE]

hip arthroplasty.

₹ Show: 20 Send to Text Sort Display Abstract

> Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer

patients had an excellent or good outcome, according to the Harris rating system, at the latest review. The two patients with poor results had obvious alternative causes for their

continuing symptoms. There was no radiological evidence of bone or prosthesis failure during the period of follow-up. CONCLUSIONS: The preliminary results are comparable with those of other authors who have examined the early results of metal-on-metal total

Tutorial New/Noteworthy

E-Utilities

LinkOut Cubby

**PubMed Services** 

Journals Database MeSH Database Single Citation Matcher

Batch Citation Matcher Clinical Queries

Related Resources Order Documents

NLM Gateway

Clinical Alerts

ClinicalTrials.gov

PubMed Central

Privacy Policy

TOXNET
Consumer Health







	VCDI	1		A.A.		of Medicine	
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	OMIM PMC	Journals Books
Search Pu	ıbMed	for		- Company of the Comp		Go C	ear
Scalch I		Limits	Previe	ew/Index	History	Clipboard	Details
About Entrez		<b>Display</b> Abstr	act	Sho	w: 20 Sort	্ ∯Send to	Text
Text Version					•		
Entrez PubM Overview Heln LEAO	led	1: Rev Chir O	EZMED	ratrice Appa	ar Mot. 1996;82(1	.):85-89.	Related Articles, Link

[15 years survival of the Mac-Kee Farrar metal hip prosthesis. Apropos of 58 cases and 4 explanted cups]

[Article in French]

Ray A.

Clinique Orthopedique du Parc de Lyon.

PURPOSE OF THE STUDY: Despite a high percentage of loosening (femoral, iliac or both), many surgeons have been surprised by some excellents results of Mac Kee Farrar prosthesis after 15, or even 20 years follow-up. MATERIALS AND METHODS: 37 patients (on 58) were reviewed with a follow-up of more than 15 years (48 hips). Among 17 cases followed for more than 20 years, with very good results (clinical and radiological), only one femoral loosening was observed. A part, 4 paired explanted implants (loosening at 18 years for 3 and at 21 years for one) were examined for a dimensional and metallurgic study. RESULTS: The results showed: no wear, very good bearing surface statement and sphericity, We never observed agressive granulomatous lesions with metallic particles (metallosis), nor wear concerning the cup. CONCLUSION: The peripheric design appears able to give a very good pressure repartition from cup to bone, allowing a homogeneous coat of cement with an equal thickness, and avoids loosening. Finally, we think that the progressive polar cavity in the cup, could have a great importance on lubrification, as an hydrokinetic reserve and micropump for synovial fluid.


PMID: 8762995 [PubMed - indexed for MEDLINE]

	·	
Display Abstract	Show: 20 Sort Send t	o Text

Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer







			₹				
Entrez	PubMed	Nucleatide	Protein Gen	ome Structure	OMIM PMC	Journals Boo	кs
Search	PubMed	for			Go	Clear	
Search		Limits	Preview/Inde	x History	Clipboard	d Details	
About Entr	ez	<b>Display</b> Abstra	ect	Show: 20 S	ort Ser	nd to Text	

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

1: Proc Inst Mech Eng [H]. 2001;215(2):153-60.

Related Articles, Links

In vitro comparison of the two hard-hard articulations for total hip replacements.

Rieker C, Konrad R, Schon R.

Sulzer Orthopedics Limited, PO Box 65, CH-8404 Winterthur, Switzerland.

Polyethylene particle disease is one of the major causes of late aseptic loosening of total hip replacement. Two hard-hard articulations (alumina-on-alumina and metal-on-metal) have been developed in Europe as an alternative to the ultra-high molecular weight polyethylene (UHMWPE) articulations. Even though these hard-hard articulations are on the market and numerous reports have been published about them, only a very limited number of studies allowing a direct in vitro comparison of the two articulations have been published so far. This paper compares in vitro these two types of articulation (alumina-on-alumina and metal-on-metal), which have been tested with a hip simulator for their tribological behaviour using exactly the same experimental methodology. This comparison shows that these two types of hard-hard articulation have very similar abrasive wear behaviour with four main features: 1. A running-in wear period (1  $\times$  10(6) cycles) gives a cumulative wear of about 20 microns with head diameters of 28 mm. 2. After the running-in wear, there is a stabilization of the linear wear behaviour with a low linear wear rate/10(6) cycles for both types of articulation. 3. The volumetric wear rate of both articulations (< 2.0 mm3/year for head diameters of 28 mm) is significantly lower than that observed for metal-on-polyethylene or ceramic-on-polyethylene articulations having the same head diameter. 4. Abrasive wear is readily apparent (indicating a mixed lubrication regime) with both types of articulation. The extremely low wear performance of these articulations is confirmed and they constitute a low-wear alternative to the UHMWPE articulations currently used.

PMID: 11382074 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort Send to Text

Write to the Help Desk

NCBI | NLM | NIH

Department of Health & Human Services
Freedom of Information Act | Disclaimer







	Pub	MILE	e <b>G</b>		ibrary Medicine		
Nucleotide	Protein	Genome	Structure	OMIMO	PMC	Journals	Books
for				kada salahan salahan salahan salahan da salah	Go] CI	ear	
Limits	Prev	iew/Index	History	(	Clipboard	De	tails
D <b>isplay</b> Abs	tract	Sho	w: 20 🕝 Sc	ort S	Send to	Text	
	•						
ELSE Sursess Retrie	VIER Marces val study	of uncen	nented me	tal-me	tal hip p		
			•				
Departm Technolo	ent of Micro ogy, Floraga	o-Technique a sse 7, A 1040	and Precision E O Vienna, Aust	ingineerir ria. reinis	ng, Vienna i ch@jic.at	University o	of
manufac 32 mont (range, 2 modern coordina	A tribologic assessment was performed on 22 metal-metal hip prostheses from a single manufacturer, following removal for early aseptic loosening after a mean service life of 32 months (range, 12-59 months). The mean linear wear rate was 7.6 microm/year (range, 2.9-12.8 microm/year). This was below the rates previously observed in other modern metal-metal combinations. A novel contour analysis technique using a coordinate measuring machine showed the mean volumetric wear rate to be 2.02 mm						
	Display Abs  1: Biomater E L S E TILLER  Retrier reviser  Reiniscl  Departm Technolo  A tribolo manufac 32 mont (range, modern	Timits Prevalent for Limits Figure 1212 Retrieval study revised for ear Reinisch G, Judman Department of Micro Technology, Floraga A tribologic assessment for Limits (range, 129-12.8 min modern metal-m	for Limits Preview/Index  Display Abstract Sho  1: Biomaterials. 2003 Mar;24(6):108  ELSEVIER  FULL SEVIER  F	Display Abstract Show: 20 Score Show	Timits Preview/Index History  1: Biomaterials. 2003 Mar;24(6):1081-91.  ELSEVIER  INTEGRATION  Reinisch G, Judmann KP, Lhotka C, Lintner F, Zwe  Department of Micro-Technique and Precision Engineerin Technology, Floragasse 7, A 1040 Vienna, Austria. reinis  A tribologic assessment was performed on 22 metal-met manufacturer, following removal for early aseptic loosen 32 months (range, 12-59 months). The mean linear weat (range, 2.9-12.8 microm/year). This was below the rate modern metal-metal combinations. A novel contour analycoordinate measuring machine showed the mean volum	of Medicine    Protein   Genome   Structure   OMIM   PMC	of Medicine  Nucleotide  Protein  Genome  Structure  OMIM  PMC  Journals  Go  Clear  Limits  Preview/Index  History  Clipboard  Def  Display  Abstract  Show:  20 Sort  Send to Text  Related A  ELSEVIER  SULFFRAGE  Retrieval study of uncemented metal-metal hip prosthese revised for early loosening.  Reinisch G, Judmann KP, Lhotka C, Lintner F, Zweymuller KA.  Department of Micro-Technique and Precision Engineering, Vienna University of Technology, Floragasse 7, A 1040 Vienna, Austria. reinisch@jic.at  A tribologic assessment was performed on 22 metal-metal hip prostheses from manufacturer, following removal for early aseptic loosening after a mean servi 32 months (range, 12-59 months). The mean linear wear rate was 7.6 microm (range, 2.9-12.8 microm/year). This was below the rates previously observed modern metal-metal combinations. A novel contour analysis technique using a

Privacy Policy

PubMed Central

PMID: 12504531 [PubMed - indexed for MEDLINE]

causative factor for early loosening.


Send to ∑ Show: 20 Text Sort Display Abstract

rate of 16.9 mg/year (range, 4.6-31.4 mg/year). The mean clearance of 39.8 microm

combinations. Evidence of abrasive, adhesive, and third-body wear was found on all bearing surfaces. The tribologic assessment did not indicate manufacturing defects as a cause of early loosening. Equally, third-body wear was too low to be considered a

(range, 30-50 microm) was within the optimal range for hard-hard bearing

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer



PMID: 7595934 [PubMed - indexed for MEDLINE]

performance for most of the categories that were examined.

 $\mathbf{z}$ Send to Text Sort Display Abstract

the stainless steel. It was found that the composite prosthesis gave the best

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer







	V			•					
Entrez	PubMed	Nucleotice	Protein	Genome	Structure	OMIM	PMC	Journals	Book
ı	PubMed	for					Go C	ear	
Search l		Limits	Preview	/index	History	С	lipboard	Det	ails
About Entre	ez	<b>Display</b> Abstra	ct	Shov	, 20 😸 So	t Š	Send to	Text	M

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ 1: Acta Orthop Scand. 2003 Aug;74(4):380-8.

Related Articles, Links

## Alternative materials to improve total hip replacement tribology.

Santavirta S, Bohler M, Harris WH, Konttinen YT, Lappalainen R, Muratoglu O, Rieker C, Salzer M.

Department of Orthopaedics and Traumatology, Helsinki University Central Hospital, Helsinki, Finland. seppo.santavirta@hus.fi

An improvement in tribology of bearing surfaces is an effective means of increasing the longevity of total hip replacement (THR). Currently, 3 approaches are available to achieve this aim: first, use of highly cross-linked UHMWPE; second, aluminum oxide ceramic bearings, and third, metal-on-metal bearings. Cross-linking reduces the wear resistance of UHMWPE markedly without impairment of other significant properties of the material. Simulator studies and some clinical long-term (10-22 years) follow-up surveys suggest an almost immeasurable wear of the highly cross-linked UHMWPEbased acetabular components during an expected clinical life span. Bioinert alumina ceramic (aluminum oxide) was introduced 3 decades ago for THR-bearing surfaces to improve performance and longevity. Alumina ceramic is entirely biostable and bioinert and has good mechanical properties. For correctly positioned alumina-on-alumina bearings, the annual linear wear rate has been reported to be 3.9 microm. Alumina heads have been successfully used in combination with polyethylene sockets, but as regards wear, the best results have been obtained with alumina-on-alumina bearings. In ceramic THR bearings, precise manufacture and contact surface geometry, including optimal clearance, are most important. For the currently available products, the component fracture risk is almost nonexistent (less than 1 per 1000). Metal-on-metal bearings were used in the early stage of THR surgery, although not all old designs were successful. More recent analyses of the early series have shown the advantages of metal-on-metal to be better and have led to a renaissance of this articulation. Initially, stainless steel was used because it was easy to manufacture and polish. Current metalon-metal bearings are based on cobalt-chromium-molybdenum alloys with varying carbon contents. Such bearings are self-polishing. Linear wear rates remain at the level of a few microm a year. An improvement in technology has increased the life span of the above three THR-bearing systems. Although the technical solutions differ considerably, they all seem to improve clearly the tribology and longevity of the THR. Each of these bearing concepts will probably permit the use of larger head sizes, to reduce the risk of impingement and luxations.

#### Publication Types:

- Review
- Review, Tutorial

PMID: 14521286 [PubMed - indexed for MEDLINE]

					•	
S NCBI	I	Pub	$\mathbb{Q}_{\epsilon}$	ed	National Library of Medicine	.M
Entrez PubMed	Nucleotide	Protein	Genome	Structure	OMIM PMC	Journals Books
Search PubMed	for	AND THE PROPERTY OF			Go Cle	<b>a</b> (
	Limits	Previe	ew/Index	History	Clipboard	Details
About Entrez	Display Abstr	act	Sho	w: 20 😽 Sort	Send to	Text
Text Version						
Entrez PubMed Overview Help   FAQ Tutorial New/Noteworthy	Ion rele	rScience ase in s	table hip	arthroplasti	es using meta	Related Articles, Links
E-Utilities		_		comparison c	etween short	anu
PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby	Savarino Baldini N Laboratori	<b>, Giunti A</b> . o di Fisiopa	n <b>i D, Ciape</b> i atologia deg	li Impianti Ortope	<b>reco M, Rotini R,</b> edici, Istituti Ortope	
Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central  Privacy Policy	The use of polyethyle conspicuousen foun rate is hig surface fir milliliter) in term (sub follow-up: implant. Consultant analyzed. The reference and Crollewith the substant pressignificant apparent population "internal" order to consists consultant medium order to	f metallic hane (PE) we us release do in the shading. Our name of the shading of the shading in the range was rels in the range was advantage advantage has, there is exposure to the shading long later Res 6	ear in total had metal ior metal ior ort-term, but the contrary, repurpose whats with metal control of metal or medium-term values (subject of metal-or a major control of metal-or control of control	ating with metalling replacement ( s from new mode of it is yet unclear it becomes neglig as to compare the cal-on-metal stable ow-up: 24 mo) and etermine whether (Co), molybdenum all patients were used to a population of (subgroup A(2)) and values. Otherwise rison with referent-metal coupling, incern about the expression of ion disseminity of metals. Coping, 2003	ible, because of the serum ion values of prosthesis (Groupd medium-term (so the ion release den (Mo) and aluminated for comparison of 27 healthy subject were significantly hand all concerce values. In conclesspecially in youngextent and duration	tal bearings. A tal bearings has um-term corrosion e continuous (nanograms per p A), in the short-ubgroup A(2); mean creased with time of um (Al) were (Group B). The ects (Group C). Co sed in comparison igher (p < 0.001) intrations were not lusion, despite the per patient of the relevant efully monitored, in ueently, to identify

Send to Text

Display Abstract

Show: 20 | Sort

PubMed Services

Journals Database MeSH Database

Clinical Queries LinkOut

Cubby

Single Citation Matcher Batch Citation Matcher

Related Resources

Order Documents

Consumer Health

PubMed Central

Privacy Policy

NLM Gateway TOXNET

Clinical Alerts ClinicalTrials.gov

S.	NCBI		Pub	Me	ed	Lib	onal orary icine N	Ĭ.M	
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	MIMO	PMC	Journals	Books
Search Pu	bMed	for			History	Clir	Go CI	<b>ear</b> De	tails
About Entrez		Limits Display Absi		ew/Index	20 <b>3</b> Son		Send to	Text	
Text Version  Entrez PubM	ed	$\square_1$ : J Biomed	Mater Res.	2002;63(5):	467-74.			Related A	urticles. Links
Overview Help   FAQ Tutorial New/Noteworth E-Utilities		Ion release in patients with metal-on-metal hip bearings i joint replacement: a comparison with metal-on-polyethyle bearings.							in total ene

Savarino L, Granchi D, Ciapetti G, Cenni E, Nardi Pantoli A, Rotini R, Veronesi CA, Baldini N, Giunti A.

Laboratorio di Fisiopatologia degli Impianti Ortopedici, Istituti Ortopedici Rizzoli, Bologna, Italy. lucia.savarino@ior.it

Polyethylene (PE) wear has been shown to be a problem in long-term joint replacement using metal-on-PE bearing. The use of metallic heads articulating with metallic cups could solve this problem: success will be enhanced if wear and corrosion of the articulating surfaces are maintained at a low level. New models with metal-on-metal bearing have been proposed, to be used mainly for young subjects: such coupling seems to have a reduced release, but it is unclear yet if the medium-term corrosion rate is really negligible or, on the contrary, it is significantly higher than in the metal-on-PE bearing. Aim of our study was the comparison of ion release in the serum of two groups of patients who had the same type of stable cementless prosthesis, but different bearing: twenty-six patients with metal-on-metal (Group A) and fifteen patients with metal-on-PE bearing (Group B) were examined. The follow-up was 14-38 months for group A and 18-34 months for group B. The serum concentration of chromium (Cr), cobalt (Co) and molybdenum (Mo) was measured. Twenty-two patients before surgery were used for comparison (Group C). The reference values were obtained from a population of twenty-two healthy subjects (Group D). Our findings indicate that metalon-metal bearings produce a significantly higher systemic release of cobalt and chromium (ng/ml) when compared with levels found in metal-on-PE, pre-surgery and reference groups. Such a high release should induce to improve the bearing materials or, at least, to study the biologic fate of metal ions and consequently their long-term effects. In such a way a risk-to-benefit ratio for the patient could be established. Copyright 2002 Wiley Periodicals, Inc. J Biomed Mater Res (Appl Biomater) 63: 467-474, 2002

PMID: 12209889 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort Send to Text
------------------	----------------------------

E-Utilities

LinkOut Cubby

TOXNET

PubMed Services Journals Database

Batch Citation Matcher Clinical Queries

Related Resources

Order Documents **NLM Gateway** 


Consumer Health

PubMed Central

Privacy Policy

Clinical Alerts ClinicalTrials.gov

MeSH Database Single Citation Matcher



replacement.

Schaffer AW, Pilger A, Engelhardt C, Zweymueller K, Ruediger HW.

Department of Occupational Medicine, University of Vienna, Austria. Andreas.Schaffer@akh-wien.ac.at

OBJECTIVE: To determine metal concentrations in blood and urine of patients who received cobalt-chromium-alloy metal on metal hip implants. METHODS: Cobalt and chromium were determined in blood and urine of 76 patients and 26 controls by electrothermal atomic absorption spectroscopy. RESULTS: A significant postoperative elevation of the metal concentrations was observed for total hip replacement patients in contrast to the control group. Twenty-nine patients exceeded the EKA (Expositionaquivalente fur Krebserzeugende Arbeitsstoffe) threshold limits for cobalt in blood and for cobalt and chromium in urine. We obtained a significant correlation between cobalt in blood and cobalt in urine (r = 0.79; p < 0.005), chromium in blood and chromium in urine (r = 0.79; p < 0.005), cobalt in blood and chromium in blood (r= 0.69; p = 0.008), and cobalt in urine and chromium in urine (r = 0.95; p = 0.004). CONCLUSION: Our findings suggest that in total hip replacements using metal-metal pairings, metal ions of the alloys are released. This release may lead to significantly elevated metal concentrations in biological fluids. Long-term studies are needed to determine the risk of metal-metal implants as a potential cause of cobalt and chromium toxicity.

Publication Types:

Clinical Trial

PMID: 10630267 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort Send to Text

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer

New/Noteworthy

**PubMed Services** Journals Database

MeSH Database Single Citation Matcher Batch Citation Matcher

Clinical Queries LinkOut

Related Resources

Order Documents **NLM Gateway** 

Consumer Health

PubMed Central

Privacy Policy

Clinical Alerts ClinicalTrials.gov

Cubby

TOXNET


E-Utilities

Details

Related Articles, Links

Text

Books



Factors correlating with long term survival of McKee-Farrar total hip prostheses.

Schmalzried TP, Szuszczewicz ES, Akizuki KH, Petersen TD, Amstutz HC.

Joint Replacement Institute at Orthopaedic Hospital, Los Angeles, CA, USA.

Clinical and radiographic evaluations of 15 McKee-Farrar hip replacements in 13 patients with followup of 21 to 26 years were performed. The average Harris hip score was 86 with no patients having a poor result. These patients outscored the age matched controls in all categories of the SF-36 health survey. All patients were community ambulators with qualitative activity levels exceeding the average for their age. Quantitative activity assessment with a pedometer in 3 patients indicated a current average of approximately 900,000 cycles per year. This represents more than 21 million cycles when extrapolated during the life of the implants. None of the femoral components were radiographically loose. One acetabular component may be loose. Osteolysis developed in 3 apparently well fixed femurs and in 1 acetabulum. There were several features of these cases that may have contributed to the long survival: (1) relatively small stature of the patients who averaged 160.5 cm (5 feet 5 inches) in height and 66.9 kg (147 lbs) in weight; (2) favorable biomechanics of the reconstruction with the hip center of rotation being medialized by an average of 6.4 mm and the femoral offset increased by an average of 4.9 mm; (3) decreased potential for neck socket impingement with an average lateral acetabular opening of 54 degrees and all components were anteverted; (4) radiolucent cement in 13 of 15 hips; and (5) no radiographically measurable wear. Previous analyses and comparisons of the clinical performance of the McKee-Farrar implant have focused on the metal on metal bearing. As has been recognized with the many variations of total hip replacement using metal on plastic hearings, there are a myriad of variables that contribute to clinical outcome. The results of this study suggest that patient selection and technical factors may contribute to the long term survival, and conversely to the failure, of McKee-Farrar implants.

Publication Types:

- Review
- · Review, Tutorial

PMID: 8769322 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort	Send to Text
Applied The Paris of the Paris	Show: I	

Tutorial

New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database

Clinical Queries LinkOut

**NLM Gateway** 

Consumer Health Clinical Alerts

ClinicalTrials.gov

PubMed Central


Privacy Policy

TOXNET

Cubby

Single Citation Matcher Batch Citation Matcher

Related Resources Order Documents



Metal on metal surface replacement of the hip. Technique, fixation, and early results.

Schmalzried TP, Fowble VA, Ure KJ, Amstutz HC.

Joint Replacement Institute at Orthopaedic Hospital, Los Angeles, CA, USA.

High volumetric wear of polyethylene plays a central role in periprosthetic bone resorption and the failure of metal on polyethylene total hip resurfacing prostheses. An assessment of technique, initial fixation, and the early results of 21 hips in 19 patients implanted with a metal on metal bearing total hip resurfacing prosthesis, 4 all cementless Wagner prostheses and 17 all cemented McMinn prostheses, is presented. Pain relief was equal to conventional total hip replacement with a better functional result with an average followup of 16 months (range, 10-25 months). The femoral component position and fixation is satisfactory in all 21 hips and there were no femoral neck notches or fractures. All 4 cementless Wagner acetabular components appear to be osseointegrated with stable interfaces. The cemented McMinn acetabular components, however, have shown progressive cement bone interface radiolucencies in 12 hips. This preliminary experience underscores the importance of obtaining secure initial fixation. There have been no problems directly attributable to the metal on metal bearing but the authors will continue to follow these hips and evaluate their performance. The metal on metal hip surface replacement procedure is in evolution. This ongoing experience will help quide total hip surface replacement component design and implantation techniques.

Publication Types:

Case Reports

PMID: 8769328 [PubMed - indexed for MEDLINE]

Display Abstract Show: 20 Sort Send to Text
---

Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer

Apr 19 2004 06:53:45

341

Help | FAQ

Tutorial New/Noteworthy

E-Utilities

PubMed Services Journals Database MeSH Database

Clinical Queries LinkOut

Cubby

Single Citation Matcher Batch Citation Matcher

Related Resources

Order Documents

**NLM Gateway** 


Clinical Alerts

ClinicalTrials.gov

PubMed Central

Privacy Policy

TOXNET
Consumer Health



Long-duration metal-on-metal total hip arthroplasties with low wear of the articulating surfaces.

Schmalzried TP, Peters PC, Maurer BT, Bragdon CR, Harris WH.

Joint Replacement Institute at Orthopaedic Hospital, Los Angeles, CA 90007, USA.

The 20-year performance of metal-on-metal hip articulations has not been reported. Five McKee-Farrar total hip prostheses and one Sivash prosthesis were obtained at revision surgery after a mean implantation time of 21.3 years. A radiographic, histologic, implant, and wear analysis was performed on these total hip implants with cobaltchrome metal-on-metal articulations. All cases were associated with femoral component loosening, but the bearing surfaces performed remarkably well. The worst case estimate of combined femoral and acetabular linear wear was 4.2 microns per year, about 25 times less than that typically seen with polyethylene. Metal particles and foreign-body inflammation were seen in all cases, but the volume of reactive tissue was small compared with what is generally seen at revision of hips with a polyethylene acetabular bearing. This may be due to a reduced particle burden or a decreased inflammatory reaction to particulate metal, or both. In addition to articular wear, other sources of metal particles included femoral neck impingement on the acetabular rim, stem burnishing, and corrosion. Prosthetic hip reconstructions can fail for many reasons, including suboptimal femoral stem and/or acetabular cup design and/or fixation. By today's standards, the McKee-Farrar and Sivash stem and acetabular component designs are suboptimal; however, after more than 20 years of use, the metal-on-metal bearing surfaces in these cases demonstrated low wear and do not appear to be the cause of failure. Recent advances in total hip arthroplasty, which include improved implant design, materials, manufacturing, and fixation, combined with a better understanding of the mechanisms of implant loosening and failure, suggest that the cobalt-chrome metal-on-metal bearing be reexamined as an alternative to polyethylene when exceptional durability is required.

Publication Types:

Case Reports

PMID: 8713913 [PubMed - indexed for MEDLINE]

Display Abstract	
------------------	--

Order Documents

Consumer Health Clinical Alerts

ClinicalTrials.gov

PubMed Central


Privacy Policy

**NLM Gateway** 

TOXNET

Details

Books




of total hip replacements. It has stimulated the development of alternative joint surfaces such as metal-on-metal or ceramic-on-ceramic implants. Since 1988 the second-generation metal-on-metal implant Metasul has been used in over 60 000 hips. Analysis of 118 retrieved specimens of the head or cup showed rates of wear of approximately 25 microm for the whole articulation per year in the first year, decreasing to about 5 microm per year after the third. Metal surfaces have a 'self-polishing' capacity. Scratches are worn out by further joint movement. Volumetric wear was decreased some 60-fold compared with that of metal-on-polyethylene implants, suggesting that second-generation metal-on-metal prostheses may considerably reduce osteolysis.

PMID: 10068001 [PubMed - indexed for MEDLINE]

Show: 20 Display Sort Send to Text Abstract

> Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer



Send to

Text

Display

Abstract

Sort

8 NCBI		Pub	Щe	d	National Library of Medicine	M.M.
Entrez PubMed	Nucleotide	Protein	Genome	Structure	OMIM PMC	Journals Books
PubMed	for	<u> </u>			(Go Cle	ear.
Search L' domed	Limits	Previe	ew/Index	History	Clipboard	Details
About Entrez	Display Abs	tract	Sho	w: 20 🕃 Sort	Send to	Text
Text Version						
	☐ 1: ) Biomed	Mater Res	2004 Jan 15	:68B(1):1-14.		Related Articles, Links
Entrez PubMed Overview Help   FAQ	inte	erScience	*			
Tutorial New/Noteworthy	Wear e	valuatio	n of coba	lt-chromium	alloy for use	in a metal-on-
E-Utilities		hip prost				
PubMed Services Journals Database MeSH Database			ckas LD, Po	oggie RA.		
Single Citation Matcher Batch Citation Matcher	Departme	ent of Ortho	paedic Surge	ery and Rehabilita	ation, The Universi	ity of Mississippi
Clinical Queries	Medical C	Center, 2500	North State	Street, Jackson,	Mississippi 39216,	, USA.
LinkOut Cubby	kstjohn@	sod.umsme	d.edu			
Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central  Privacy Policy	early dev research by modification on-metal systems componed data public of a metal disk testion wear research metal hij was investigation.	continues to continues to continues to continues to continues to continue the continue that has been allic femoral lar metal. The continue for using showed istant, and to codeled the lets of metalmetric wear 180 times to from either rized and continue prostheses estigated hole criodicals, In continue to continue the continue the continue to continue the conti	which has been attempt to olymer procesuples for hip e couple of reperformed mead compositions study involved and the allocations alloy was conspolyethylof the metal that for the reperformed with a semanufactured sufficient of a sufficient of sufficient of the sufficient of suffici	reduce the wear reduce the wear ressing, there is a prostheses. Wean etal or ceramic land reported, but hip-simulator we conent with an acception of the hip-bet with the higher of the samples to make the samples to make the samples to make the particles retrievabler. Based upon the higher commise to justify	of polyethylene journement interest in testing of total has been polymeriat, until recently, the ar studies involving trabular cup compositro wear resistance of primarily in the earing couple. The (0.25%) carbon confor testing in a has been polyethylene testing or hip simple. Polyethylene testing or hip simple from periprostit the results of this carbon cobalt/chi	c acetabular there has been little ig the combination osed of the same or ce of two c carbon content, as results of pin-on- content was more hip-simulator system, lly. Comparison of inples showed that fer 5,000,000 cycles and metal particles ulator testing were netic tissues by other

Send to Text

Display Abstract

Show: 20 Sort



Privacy Policy





Entrez	PubMed	Nucleotide	Protein	Genome	Structure	OMIM	PMC	Journals	B <sub>'</sub>	
Search	PubMed	<b>⊡</b> for	A LEGISTRY I PROMING I PRILITY Y PROST	gg/frameg safe/ I wildow N cost if all MISS/2 h. I Meditioners N o N 140	en den in , min, land 100m, aver sammer det sessionalistich, but en preparet zu im	cooperation to the property of the fields of	Go	Clear		
**		Limits	Previ	ew/Index	History	Cli	pboard	Def	tails	
About Enti	ez	Display	act	Sho	ow: 20 🕞 S	ort §	<b>√</b> Send	to Text	je K	
Text Versi	on	□1: J Arthro	olasty.	1998 Apr;	13(3):343-5	; discus	sion	Related A	rticles	
Entrez Proverview Help   FACT Tutorial New/Note E-Utilities	2	345-6.  Commer  • <u>J A</u>	nt on: rthropl	asty. 1998	Apr;13(3):3	339-43;	discuss			
D l- N 1 l	Caminon						mp ur o	opano	-5	
PubMed Journals I MeSH Da Single Cit Matcher	Database tabase ation	Streiche		ising mei	tal-on-me	.a				
Batch Cita Matcher	ation	Stryker, Adliswil, Switzerland.								
Clinical Q LinkOut Cubby	ueries	Publication Types:  • Comment								
Related Resource Order Doo		PMID: 9	959064′	7 [PubMed	l - indexed f	for MEI	OLINE]			
NLM Gate TOXNET Consume Clinical A ClinicalTr PubMed (	r Health lerts ials.gov	Display Abstr	act	Sh	ow: 20 🔁 S	Sort	Send	to Text		

Write to the Help Desk

NCBI | NLM | NIH

Department of Health & Human Services

Freedom of Information Act | Disclaimer

Entrez PubMed Overview Help | FAQ

PubMed Services Journals Database MeSH Database

Single Citation Matcher Batch Citation Matcher Clinical Queries

Related Resources Order Documents

**NLM Gateway** 

Consumer Health Clinical Alerts

ClinicalTrials.gov

PubMed Central

Privacy Policy

TOXNET

Tutorial New/Noteworthy


E-Utilities

LinkOut

Cubby







			7 ~	₹			
Entrez	PubMed	Nucleotide	Protein Genon	ne Structure	OMIM PMC	Journals 8	ooks
Search Pu	bMed	for			Go C	lear	
Jean Cit /		Limits	Preview/Index	History	Clipboard	Details	
About Entrez		Display Abstrac	ct 👼	Show: 20 So	nt ⊋ Sendit	o Text	
Text Version							

☐ **1:** J Arthroplasty. 1997 Oct;12(7):819-24.

Related Articles, Links

## Progressive bilateral pelvic osteolysis in a patient with McKee-Farrar metal-metal total hip prostheses.

Szuszczewicz ES, Schmalzried TP, Petersen TD.

Harbor-UCLA Medical Center, Torrance, California, USA.

As accumulating evidence indicates that polyethylene plays a central role in periprosthetic osteolysis, there is a renewed interest in alternatives such as metal-metal bearings. Several long-term studies report encouraging results with the McKee-Farrar total hip arthroplasty, but there is a paucity of data on the incidence, severity, and pathogenesis of osteolysis in metal-metal bearing total hip arthroplasties. This study presents a patient who had progressive bilateral pelvic osteolysis associated with his McKee-Farrar metal-metal total hip prostheses. His left hip was revised after 13.5 years of service. The tissues revealed no gross metal staining and fewer inflammatory constituents than are typically found in metal-polyethylene bearing hips. His right hip was still functioning after 22.5 years of service, although the acetabular component was loose by that point. An arthrogram of this hip demonstrated communication of the joint with the iliac osteolysis. The development of osteolysis in both hips followed a pattern similar to that seen in metal-polyethylene total hip arthroplasties. Bearing wear could not be detected in either of the hips. Accumulating evidence indicates that particulate debris of appropriate size and number is capable of fueling periprosthetic inflammation. Specific to this study, consideration should be given to particles of cobalt-chromium alloy, polymethyl methacrylate bone-cement, and barium sulfate. Other factors that should be considered are increased joint fluid pressure, soluble inflammatory mediators, and the effective joint space. When bone becomes part of the effective joint space, it is exposed to particulate debris, soluble factors, and potentially increased joint fluid pressures, which may promote localized bone resorption. It must be kept in mind that the development of osteolysis is multifactorial. Although bearings with better wear characteristics are desirable, the elimination of polyethylene will not eliminate osteolysis.

Publication Types:

· Case Reports

PMID: 9355013 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort	Send to Text
------------------	---------------	--------------







······			•	• •					
Entrez	PubMed	Nucleotide	Protein C	Senome	Structure	OMIMO	PMC	Journals.	Book
Search P	ubMed	for					Go C	lear	
Jear Cil Lune		Limits	Preview/Ir	ndex	History	C	lipboard	De	etails
About Entrez		Display Abstra	ct	Shov	v: 20 🙀 So	rt 🔀	Send to	o≫ Text	Ä

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

**1:** Orthopade. 1997 Feb;26(2):142-51.

Related Articles, Links

|SpringerLink

[Changes in shape of the McKee-Farrar hip endoprosthesis]

[Article in German]

Tager G, Euler E, Plitz W.

Chirurgische Klinik und Poliklinik, Klinikum Innenstadt, Ludwig-Maximilians-Universitat Munchen.

The still unsolved problem of aseptic loosening in total hip arthroplasties with identification of polyethylene wear particles as one of its major causes, has led to reintroduction of metal-to-metal articulations, as indicated by a few good clinical longterm results with all-metal McKee-Farrar arthroplasties. In this paper, data on 145 patients from a population of more than 1400, all with implanted McKee-Farrars, who underwent revision surgery for aseptic loosening, are collected and analysed for dependence of duration to brands of the implants and position of the cups. The surface of each of 55 revised implants was measured using a 3-D device. The results showed no interdependence between time of loosening, brand inclination of the cup and deviation in shape of ball and cup. Additionally, the deviations in shape were slight.

PMID: 9157354 [PubMed - indexed for MEDLINE]

Display Abstract	Send to Text	<b>Y</b>
Part Control of the C	Show: 1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-	

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer

Privacy Policy

S NCBI	į	Pub	Me	d		ional brary licine	NLM	
Entrez PubMed	Nucleotide	Protein	Genome	Structure	OMIM	PMC	Journals	В
Search PubMed	☑ for l	Provi	ew/Index	History	- week commence of the commenc	Go C oboard	Clear De	etails
About Entrez	Display Absti			w: 20 So		Send t	o Text	
Text Version	☐ 1: J Bone Jo	e A de Composito d			74-80.	Rel	ated Articles	, Links
Entrez PubMed Overview Help   FAQ Tutorial New/Noteworthy E-Utilities	The ris	xt crticle co .ejbjs.org k of can	cer follow	ing total hip		e arthr	oplasty.	
PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby	Joint Rep 90007, U Publicati • Re	placement ISA. on Types view	:	Orthopaedic	Hospital,	Los Anş	geles, CA	
Related Resources Order Documents NLM Gateway TOXNET	<u> </u>	view Lite 1379749		ndexed for M	EDLINE	]		
Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central	Display Abst	ract	Sho	w: 20 - So	nt S	Send	to Text	

Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer







S NGBI		Pubywed				of Medicine				
Entrez PubMed	Nucleotide	Protein	Genome	Structure	MIMO	PMC	Journals	Books		
Search PubMed	for			· 45 (2014) (		Go C	ear			
	Limits	Previ	ew/Index	History	CI	ipboard	Det	ails		
About Entrez	and the second s					greenkrennij (nistillijk)				
	Display Abst	tract	Sho	w: 20 🛨 So	rt 🔻	Send to	Text			
Text Version										
Entrez PubMed	☐ <b>1</b> : Health Te	chnol Asses	s. 2002;6(15	):1-109.			Related A	rticles, Links		
Overview	FUL	LTEXT								
Help   FAQ Tutorial	Δ syste	matic re	view of t	he effective	eness a	nd cost	-effectiv	eness		
New/Noteworthy E-Utilities	of meta	al-on-me	tal hip re	surfacing a	arthrop	lasty fo	r treatm	ent of		
	hip dis		•	_	-	-		•		
PubMed Services Journals Database	•									
MeSH Database Single Citation Matcher	Vale L, V	Vyness L,	McCormack	K, McKenzie	L, Brazz	elli M, St	earns SC.			
Batch Citation Matcher Clinical Queries	Health Se	ervices Rese	arch Unit. In	stitute of Appli	ed Health	Sciences,	University	of		
LinkOut	Aberdeen									
Cubby		_								
Related Resources Order Documents	Publicatio • Re									
NLM Gateway		view, Acade	emic							
TOXNET Consumer Health		, , , , , , , , , , , , , , , , , , ,								
Clinical Alerts ClinicalTrials.gov	PMID: 12	:137721 [Pu	bMed - inde	xed for MEDLIN	NE]					
PubMed Central					<u></u>					
Privacy Policy				[00 Bel ]			<b>.</b>	7778		
	Display Abs	tract	Sho	w. 20 🔀 So	rt 💉	Send to	Text			


Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer

Apr 19 2004 06:53:45

350







			*	•					
Entrez	PubMed	Nucleotide	Protein G	enone	Structure	MIMO	PMC	Journals	Book
Search P	ubMed	for					Go CI	eag	
Jearen		Limits	Preview/Ind	dex	History	С	lipboard	Deta	ails
About Entrez		<b>Display</b> Abstra	ct	Show	20 👺 So	rt 🔀	Send to	Text	

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents **NLM Gateway** TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ 1: Orthopedics. 1991 Feb;14(2):137-42.

Related Articles, Links

## Cancer risk after Mckee-Farrar total hip replacement.

Visuri T, Koskenvuo M.

Department of Public Health, University of Helsinki, Finland.

Cancer incidence in 433 McKee-Farrar total hip replacement patients, operated on between 1967 and 1973, was examined for 5729 person-years, to the end of 1981. The expected number of natural deaths was slightly higher than observed, suggesting some selection of the operated patients. The risk of total cancer incidence did not increase, but the risk for site-specific cancer did because there were no cases of kidney or bladder cancer, or rare forms of cancer. The risk of leukemias and lymphomas increased, and the risk of breast cancer decreased; these results were surprisingly similar to those of a study from New Zealand. This study concluded that patients with total hip prostheses have a cancer morbidity differing from the general population. The role of chrome-cobalt-molybdenum alloy in carcinogenesis requires further investigation.

PMID: 2008381 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 💆 Sort 😤	Send to Text
Live Finell 1 consider And 1		

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer

Tutorial New/Noteworthy

E-Utilities

PubMed Services Journals Database MeSH Database

Clinical Queries LinkOut

**NLM Gateway** 

Consumer Health Clinical Alerts

ClinicalTrials.gov

PubMed Central

Privacy Policy

TOXNET


Cubby

Single Citation Matcher Batch Citation Matcher

Related Resources Order Documents







	ACDI	*	anth	100	of Medicine	NLM
Entrez	PubMed	Nucleotide	Protein Genom	e Structure	OMIM PMC	Journals Books
Search Pu	bMed	for			Go S	Clear
		Limits	Preview/Index	History	Clipboard	Details
About Entrez		Display Abstrac	et 💥	Show: 20 🖫 Sor	t 🔀 Send	to Text
Text Version						

Entrez PubMed
Overview
Help | FAQ

Canada visit offer metal on metal

Related Articles, Links

# Cancer risk after metal on metal and polyethylene on metal total hip arthroplasty.

Visuri T, Pukkala E, Paavolainen P, Pulkkinen P, Riska EB.

Central Military Hospital, Helsinki, Finland.

The incidence of cancer after metal on metal total hip arthroplasty (McKee-Farrar) and polyethylene on metal total hip arthroplasty (Brunswik, Lubinus) was compared with that of the general population in Finland. The mean followup time for the patients who had metal on metal total hip arthroplasty was 15.7 (9092 person years) and for the patients who had polyethylene on metal total hip arthroplasty it was 12.5 years (19,846 person years). One hundred thirteen malignant cancers were observed in patients who had metal on metal total hip arthroplasty and 212 were observed in patients who had polyethylene on metal total hip arthroplasty. The standardized incidence ratio for all cancers of the metal on metal arthroplasty group was 0.95 (95% confidence limits 0.79-1.13) and that of the polyethylene on metal arthroplasty group was 0.76 (95% confidence limits 0.68-0.86). The risk of total cancer in the patients who had metal on metal total hip arthroplasty was 1.23-fold compared with that of the patients who had polyethylene on metal total hip arthroplasty. Both groups had significantly less lung cancer than the general population: the leukemia incidence in the patients who had metal on metal total hip arthroplasty was slightly increased (observed to experienced 7/3.03, standardized incidence ratio 0.61; 95% confidence limits 0.17-1.56). The leukemia rate of the patients who had metal on metal total hip arthroplasty was 3.77fold compared with that of the patients who had polyethylene on metal total hip arthroplasty, but this difference was not statistically significant. No sarcomas were observed at the site of the prosthesis. The incidence of the other forms of cancers did not differ significantly from those in the general population. The observed variation in the incidence of different cancers among patients who had total hip arthroplasty compared with the general population suggests that factors other than total hip arthroplasty play a major role in the origin of cancer.

Publication Types:

- Review
- Review, Tutorial

PMID: 8769342 [PubMed - indexed for MEDLINE]

Display Abstract	20 🗐	Sort &	Send to	Text	J
Display Abstract	Show: 20	Soft 📆	Selid to	Text	







			,				
Entrez	PubMed	Nucleotide	Protein Gen	ome Structure	OMIM PI	MC Journals	Books
Search	PubMed	for			] Go	Clear	
Deare		Limits	Preview/Inde	ex History	Clipboa	ard Det	ails
About Entre	ez	Display Abstra	not 5	20 😾 S	Sort 🚽 S	end to Text	
<b>-</b>	<del></del> -	Display Abstra	ict	Show: Leville	ou 🧣 🖂	CONTRACTOR OF THE PROPERTY OF	<u> </u>

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents NLM Gateway TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

☐ **1:** Clin Orthop. 2000 Oct;(379):123-33.

Related Articles, Links



# Medium-term results of a modern metal-on-metal system in total hip replacement.

Wagner M, Wagner H.

Orthopaedic Department, Zeisigwaldkliniken Bethanien, Chemnitz, Germany.


Since 1988, metal-on-metal articulation from cobalt-chromium-molybdenum alloy was reintroduced into hip arthroplasty as an alternative to metal-on-polyethylene or ceramicon-polyethylene components. Modular joint surfaces were developed for the second generation metal-on-metal articulation using newly introduced and proven prosthetic implants. Since 1990, 78 patients with 78 uncemented total hip replacements were followed up in a prospective study. The mean followup was 60 months. Three patients were lost to followup. The average age of the patients at the time of surgery was 48.8 years. Thirty-three patients had been operated on previously. No early infections occurred; one late infection occurred after 3 years. Dislocation of the prosthesis occurred in one patient who was lost to followup. In two patients ectopic ossifications were removed 17 and 27 months postoperatively. At revision surgery no metallosis could be identified. At the last followup examination, the Harris hip score was 96.8 points on average. There was no evidence that the metal-on-metal articulation gave rise to new problems or complications. Metal-on-metal articulation reduced wear considerably in the authors' previous experience. It is hoped that foreign body reactions are reduced significantly so that an alternative for total hip replacement in younger and active patients will be available.

PMID: 11039799 [PubMed - indexed for MEDLINE]

	and the second s	
Display Abstract	Show: 20 Sort Send to Text	

Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer

Related Articles, Links



Entrez PubMed Overview

Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources
Order Documents
NLM Gateway
TOXNET
Consumer Health
Clinical Alerts
ClinicalTrials.gov
PubMed Central

Privacy Policy

Preliminary results of uncemented metal on metal stemmed and resurfacing hip replacement arthroplasty.

Wagner M, Wagner H.

1: Clin Orthop. 1996 Aug;(329 Suppl):S78-88.

Orthopaedic Hospital, Wichernhaus, Schwarzenbruck, Nuremburg, Germany.

Seventy uncemented stemmed total hip replacements and 35 uncemented surface replacements with all metal Metasul articulating surfaces were followed up in a prospective study. There was no evidence that this metal on metal articulation causes new problems or complications that were not known already from other polyethylene-aluminum oxide ceramic articulating combinations. The results of 64 of 70 patients could be assessed as excellent and good. When tissue samples obtained during 2 reoperations for ectopic ossification were examined histologically, there was no light microscopic evidence of metal particles. In these cases, aseptic loosening seemed to be due to the lack of initial fixation with the original femoral component design, and was not related to the use of the Metasul bearing. The metal on metal articulation reduces the production of particles considerably according to experience to date. It is therefore hoped that foreign body reactions due to wear particles will be significantly reduced. The results support the continued investigation of metal on metal joint replacements for younger, active patients.

PMID: 8769325 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort Send to Text

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer

Apr 19 2004 06:53:45

354



ClinicalTrials.gov PubMed Central

Privacy Policy






9	NGBL		Pub	Me	u	of Medicine	NLM
Entrez	PubMed	Nucleotide -	Protein	Genome	Structure	OMIM PMC	Journals Books
Search F	PubMed	for				Go	Clear
Jenicii j		Limits	Prev	view/Index	History	Clipboard	Details
About Entrez	!	Display Abs	tract	Show	w: 20 Sort	Send	Ito Text 🔀
Text Version							
Entrez Pub Overview	Med	🗆 <b>1:</b> Orthopad	e. 1989 Se	p;18(5):370-6	• .		Related Articles, Link
Help   FAQ		[Polve	thvlene	wear and	late loosenir	ng of a tota	l prosthesis of
Tutorial New/Notewo	rthy	the hir	ioint. N	lew persp	ectives for m	etal/metal	pairing of the
E-Utilities			e and he				
PubMed Se		•		-			
Journals Dat MeSH Datat Single Citati	oase	[Article ii	n German]				
Batch Citation Clinical Que	on Matcher	Weber I	BG, Fiecht	er T.			
LinkOut Cubby		Twenty-r	nine years'	experience wi	th total hip repla	cement surger	y, including continuing
		quality o	ontrol and i	follow-up stud	lies, indicate that	: late loosening	after 10-15 years is
Related Re Order Docui		related t	o wear of t	he polyethyler	ne cup. Metal-to-	metal prostnes	es implanted 20 years sening. In conclusion,
NLM Gatew TOXNET	ay	ago ano polvethy	longer, on lene does r	not last long e	nough when com	pared to bioco	mpatible precision-
Consumer H		made ar	tificial meta	al-metal hip jo	ints.	1	•
CHILICAL MICH	(J						

PMID: 2812770 [PubMed - indexed for MEDLINE]

Show: 20 Send to Text Sort Display Abstract

> Write to the Help Desk
> NCBI | NLM | NIH
> Department of Health & Human Services Freedom of Information Act | Disclaimer



Order Documents


Consumer Health Clinical Alerts ClinicalTrials.gov

PubMed Central

Privacy Policy

NLM Gateway TOXNET





All Andrews	NOUL	•		1		OI ME	cicine E	LAVE	
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	OMIM	PMC	Journals	Books
Search [F	PubMed	for		CHARLES AND A CONTRACTOR OF			Go C	lear	
		Limits	Prev	riew/Index	History	С	lipboard	Det	ails
About Entrez	PubMed  Limits Preview/Index History Clipboard  Display Abstract Show: 20 Sort Send to Text  Version  It version  Total hip joint replacement using a CoCrMo metal-metal pairing.  Whoteworthy Unitities  Display Abstract Show: 20 Sort Send to Text  Total hip joint replacement using a CoCrMo metal-metal pairing.  Weber BG, Semlitsch MF, Streicher RM.  Orthopadie am Rosenberg, St. Gallen/Heiden, Switzerland.  Late loosening in total hip replacement is explained by a foreign body reaction connective tissue along the bone-implant interface due to polyethylene debrise.	Text	¥						
Text Version		,							
Overview Help   FAQ Tutorial		Total h	ip joint				o meta		ticles, Link: liding
	=	Weber B	G, Semlit	sch MF, Stre	eicher RM.				
MeSH Datab	oase	Orthopad	lie am Rose	enberg, St. G	allen/Heiden, S	witzerlan	d.		
		connectiv contrast.	e tissue al metal-met	ong the bone al prostheses	-implant interfa of the McKee	ace due to type impl	o polyethy anted in th	lene debris. ne sixties, m	In ay still
Related Re	esources	work pen	rectly today	y without any	signs or osteo	iysis along	שונים שנו בי	-cement mo	=11000;

PMID: 8336059 [PubMed - indexed for MEDLINE]

Display Abstract Show: 20 Sort Send to Text

proving that metal-metal pairing is superior to prostheses with polyethylene cups and

proving also that cement anchorage may be adequate for fixation. A new metal-metal

total hip joint is presented, that has been implanted 90 times between 1988 and 1991.

Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer

Apr 19 2004 06:53:45

356







130-50 mm		Ł	7				
Entrez	PubMed	Nucleotide	Protein Ge	enome Stru	octure OMIN		Journals Book
Search P	ubMed	for				Go Cle	ar
Jearen I		Limits	Preview/Inc	dex His	story	Clipboard	Details
About Entrez		Display Abstra	nct	20	Sort	→ Send to	Text

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents **NLM Gateway** TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

1: Clin Orthop. 1996 Aug;(329 Suppl):S69-77.

Related Articles, Links

## Experience with the Metasul total hip bearing system.

Weber BG.

Orthopadie am Rosenberg, Heiden, Switzerland.

The author and Sulzer Medical Technology Ltd, Switzerland, have developed a new generation of metal on metal bearing total hip joints. The design is different than the McKee type prostheses in that the cobalt chrome alloy heads and cups (Metasul) are of the highest precision with controlled loose fitting. These allow low friction and low wear of approximately 5 micrometers per year. It is anticipated that debris related late loosening will be avoided by the use of this design. Approximately 30,000 Metasul hearings have been produced. The first 110 Weber metal on metal hip implants have been analyzed. No adverse effects from the wear of the new metal on metal components have been noted in this series.


PMID: 8769324 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort ✓	Send to Text

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer







		•	,	**		
Entrez	PubMed	Nucleotide	Protein Genor	ne Structure	OMIM PN	and became to the control of the
Search .	PubMed	for			Go	Clear
Dearent,	N	Limits	Preview/Index	History	Clipboa	ard Details
About Entre	ez	<b>Display</b> Abstra	ect 🔣	Show: 20 🔀	Sort Se	end to Text

About Entrez

Text Version

Entrez PubMed Overview Help | FAQ Tutorial New/Noteworthy E-Utilities

PubMed Services Journals Database MeSH Database Single Citation Matcher Batch Citation Matcher Clinical Queries LinkOut Cubby

Related Resources Order Documents **NLM Gateway** TOXNET Consumer Health Clinical Alerts ClinicalTrials.gov PubMed Central

Privacy Policy

1: Clin Orthop. 1996 Aug;(329 Suppl):S160-86.

Related Articles, Links

Wear behavior and histopathology of classic cemented metal on metal hip endoprostheses.

Willert HG, Buchhorn GH, Gobel D, Koster G, Schaffner S, Schenk R, Semlitsch

Department for Orthopaedics, University of Gottingen, Germany.

The authors reviewed their collection of retrieved all metal hip joints (9 McKee-Farrar, 7 Muller, and 3 Huggler type prostheses) and tissues from the joint capsules and implant beds. The amount of wear was measured, and the total volume was calculated. The tissues were analyzed by atomic absorption spectral analysis or inductively coupled plasma mass spectrometry and examined by light and scanning electron microscopy. The size of particles was measured with a texture analysis system. The articulating surfaces showed many delicate scratches which represent normal wear. The calculated annual wear averaged approximately 5 mm3 per year, which is low compared with polyethylene. The cellular reaction to metal wear particles was regarded as mild. The cellular reaction to scattered and worn bone cement was always more pronounced than to metallic debris. Scanning electron microscopy confirmed the irregular shapes and mostly submicron size of the metal particles. The analytically detected metal content of the periarticular tissue was relatively low and in accordance with the wear measurements from the articulating surfaces. The excess of chromium in the tissues is discussed in the light of the elimination of cobalt as well as the relation between elements representing either corrosion products or elements still bound in wear particles.


PMID: 8769333 [PubMed - indexed for MEDLINE]

Display Abstract	Show: 20 Sort Send to Text	4

Write to the Help Desk NCBI | NLM | NIH Department of Health & Human Services Freedom of Information Act | Disclaimer







\$55W	NCBI	i	Pub	Me	u	of Med	licine 🖪	LM	
Entrez	PubMed	Nucleotide	Protein	Genome	Structure	OMIM	PMC	Journals	Br
Search	PubMed	for	campy construction to detail of the construction of the last	manage to the second strategy to the comment date to the comment of the comment of the comment of the comment	iongropper mysses dat Messages com som en Messages co Jacobs de Message	Notice and analysis to the graphest convergency.	Go C		
•	Company of the Company Carlo State Conference on the Company of th	Limits	Previ	ew/Index	History	Cli	pboard	De	etails
About Ent	rez	Display Abst	ract	Sho	w: 20 - Sc	ort 🖹	Send to	Text	I
Text Versi	ion	□1: Clin Orth	юр. 1996	Dec;(333):4	-14.		Rela	ated Articles	, Links
Entrez P Overview Help   FAG Tutorial New/Note E-Utilities	Q	endopr	reaction ostheses HG, Sem	•	and meta	llic wear	produ	cts of joi	int
Journals I MeSH Da Single Cit	tabase ation Matcher ation Matcher	<ul><li>Bi</li><li>Cl</li><li>Hi</li></ul> Personal	on Types ography assical A storical A Name as	rticle Article			· .		
Related Order Do NLM Gate TOXNET Consume Clinical A	eway er Health	• Se	illert HG emlitsch N 3981878 [		dexed for M	EDLINE]			
ClinicalTr PubMed Privacy F	Central	Display Abs	tract	Sho	ow: 20 - S	ort	Send	to Text	Manager Annager

Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer

**NLM Gateway** 

Clinical Alerts ClinicalTrials.gov

Consumer Health

PubMed Central

Privacy Policy

TOXNET

ENGIGZ FODMICO	•				-05
S NCBI	F	op M	ed	National Library of Medicine	M
Entrez PubMed	Nucleolide	Protein Genome	Structure	OMIM PMC	Journals Books
PubMed	for			Go Cle	ac
Search I ubined	Limits	Preview/index	History	Clipboard	Details
About Entrez	Display. Abstra	act 🛣 s	how: 20 🕃 So	rt 😕 Send to	Text 🔀
Text Version					
Entrez PubMed	1: Clin Orthop.	. 1989 Sep;(246):39	-47.		Related Articles, Links
Overview Help   FAQ Tutorial New/Noteworthy	Experier	nce with total l	nip arthropla erence to lon	sty in Greece, t g-term results v	he first 20 with the
E-Utilities		arrar techniqu			
PubMed Services Journals Database MeSH Database		AL, Patikas AF.			
Single Citation Matcher Batch Citation Matcher Clinical Queries	Asklepion F	Red Cross Orthopae	dic Hospital, Athe	ens, Greece.	
LinkOut Cubby	porformed	in veny small numbe	ers However evi	ireece (1966-1967) are en after the difficult e	early period,
Related Resources Order Documents	statistics at	re low compared to	other countries.	An estimate brings th	ie total number of

Total hip arthroplasty was introduced early in Greece (1900-1907) and was intidary performed in very small numbers. However, even after the difficult early period, statistics are low compared to other countries. An estimate brings the total number of operations during a 20-year period to 9000 with a rate in recent years of 1000 per year. An early series of 143 arthroplasties (122 patients), mainly of the McKee-Farrar metal-to-metal technique, was reviewed. A final group of 52 arthroplasties, all primary prostheses of the McKee-Farrar type, were assessed with a follow-up period ranging from 12 to 20 years postoperatively. In the surviving cases, 53% were pain-free, and, in 79%, useful motion was maintained. The roentgenographic results were less satisfying but a fair roentgenographic picture did not preclude a good or very good clinical and functional outcome. Although the metal-to-metal technique now appears to be more of historic value, long-term results with this type of implant offer grounds for comparison with current cemented techniques.

PMID: 2766621 [PubMed - indexed for MEDLINE]

Display Abstract Show: 20 Sort Send to Text

Write to the Help Desk
NCBI | NLM | NIH
Department of Health & Human Services
Freedom of Information Act | Disclaimer



# DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service Food and Drug Administration Memorandum

			_
From:	Reviewer(s) - Name(s) tristaphen thack		
Subject:	510(k) Number		
То:	The Record - It is my recommendation that the subject 510(k) Notific	ation:	÷
[] _1 _1	Refused to accept.  Requires additional information (other than refuse to accept).  Its substantially equivalent to marketed devices.  NOT substantially equivalent to marketed devices.  Other (e.g., exempt by regulation, not a device, duplicate, etc.)		
	Is this device subject to Section 522 Postmarket Surveillance? Is this device subject to the Tracking Regulation? Was clinical data necessary to support the review of this 510(k)? Is this a prescription device? Was this 510(k) reviewed by a Third Party? Special 510(k)? Abbreviated 510(k)? Please fill out form on H Drive 510k/boilers	☐YES ☐YES ☐YES ☐YES ☐YES ☐YES ☐YES ☐YES	NO NO NO NO NO
	Truthful and Accurate Statement Requested Enclosed  A 510(R) summary OR A 510(k) statement  The required certification and summary for class III devices  The indication for use form		
	Combination Product Category (Please see algorithm on H drive 510k	/Boilers) <u>()</u> ([	
	Animal Tissue Source YES NO Material of Biological O	rigin 🛮 YES	M. M.
□ No	The submitter requests under 21 CFR 807.95 (doesn't apply for SEs):  Confidentiality		eding 90 d
Predic	cate Product Code with class: Additional Product Code(s) with	h panel (optiona	1):
.Kw	A, DL ClasTL  ORDB	(Date) (Date)	

## Internal Administrative Form

		VEC	NO
		YES	NO
1.	Did the firm request expedited review?		
2.	Did we grant expedited review?		<b>.</b>
3.	Have you verified that the Document is labeled Class III for GMP		
	purposes?		
4.	If, not, has POS been notified?		
5.	Is the product a device?		
6.	Is the device exempt from 510(k) by regulation or policy?		
7.	Is the device subject to review by CDRH?	<u> </u>	<b>_</b>
8.	Are you aware that this device has been the subject of a previous NSE		
·	decision?		
9.	If yes, does this new 510(k) address the NSE issue(s), (e.g.,		
	performance data)?		
10	. Are you aware of the submitter being the subject of an integrity		
	investigation?		
11	. If, yes, consult the ODE Integrity Officer.		
12	. Has the ODE Integrity Officer given permission to proceed with the		
	review? (Blue Book Memo #I91-2 and Federal Register 90N0332,		
	September 10, 1991		

## SCREENING CHECKLIST FOR ALL PREMARKET NOTIFICATION [510(k)] SUBMISSIONS

510(k) Number: <u>K04 284</u>						
The cover letter clearly identifies the type of 510(k) submission as (Check the appropriate box):						
	Special 510(k)	-	Do Sections 1 and 2			
	Abbreviated 510(k)	-	Do Sections 1, 3 and 4			
<b>9</b>	Traditional 510(k) or n	ıo identi	fication provided - Do Sections 1 and 4			

## Section 1: Required Elements for All Types of 510(k) submissions:

	Present or Adequate	Missing or Inadequate
Cover letter, containing the elements listed on page 3-2 of the	<u> </u>	
Premarket Notification [510)] Manual.		
Table of Contents.		
Truthful and Accurate Statement.		
Device's Trade Name, Device's Classification Name and Establishment Registration Number.		
Device Classification Regulation Number and Regulatory Status (Class I, Class II, Class III or Unclassified).		
Proposed Labeling including the material listed on page 3-4 of the Premarket Notification [510)] Manual.		
Statement of Indications for Use that is on a separate page in the premarket submission.		
Substantial Equivalence Comparison, including comparisons of the new device with the predicate.	V	
510(k) Summary or 510(k) Statement.		
Description of the device (or modification of the device) including diagrams, engineering drawings, photographs or service manuals.		
Identification of legally marketed predicate device. *		
Compliance with performance standards. * [See Section 514 of the Act and 21 CFR 807.87 (d).]		
Class III Certification and Summary. **	1	
Financial Certification or Disclosure Statement for 510(k) notifications with a clinical study. * [See 21 CFR 807.87 (i)]	V	
510(k) Kit Certification ***		1

<sup>\* -</sup> May not be applicable for Special 510(k)s.

\*\* - Required for Class III devices, only.

\*\*\* - See pages 3-12 and 3-13 in the Premarket Notification [510)] Manual and the Convenience Kits Interim Regulatory Guidance.

Section 2: Required Elements for a SPECIAL 510(k) submission:

	Present	Inadequate or Missing
Name and 510(k) number of the submitter's own, unmodified		
predicate device.		
A description of the modified device and a comparison to the		
sponsor's predicate device.		
A statement that the intended use(s) and indications of the	,	1
modified device, as described in its labeling are the same as the		
intended uses and indications for the submitter's unmodified		1
predicate device.		
Reviewer's confirmation that the modification has not altered the		
fundamental scientific technology of the submitter's predicate	争。 第158章	
device.		0.600,000,000,000,000
A Design Control Activities Summary that includes the following		10000
elements (a-c):	rain established	(5°5) (10°5) (1°5) (1°5)
a. Identification of Risk Analysis method(s) used to assess the		
impact of the modification on the device and its components, and		
the results of the analysis.		
b. Based on the Risk Analysis, an identification of the required	]	
verification and validation activities, including the methods or		
tests used and the acceptance criteria to be applied.		
c. A Declaration of Conformity with design controls that includes		
the following statements:		
A statement that, as required by the risk analysis, all		
verification and validation activities were performed by the		
designated individual(s) and the results of the activities		
demonstrated that the predetermined acceptance criteria were		
met. This statement is signed by the individual responsible		J
for those particular activities.		
A statement that the manufacturing facility is in conformance		
with the design control procedure requirements as specified		
in 21 CFR 820.30 and the records are available for review.		
This statement is signed by the individual responsible for		
those particular activities.		<u> </u>

Section 3: Required Elements for an ABBREVIATED 510(k)\* submission:

	Present	Inadequate or Missing
For a submission, which relies on a guidance document and/or special control(s), a summary report that describes how the guidance and/or special control(s) was used to address the risks associated with the particular device type. (If a manufacturer elects to use an alternate approach to address a particular risk,		
sufficient detail should be provided to justify that approach.		
For a submission, which relies on a recognized standard, a declaration of conformity [For a listing of the required elements of a declaration of conformity, SEE Required Elements for a Declaration of Conformity to a Recognized Standard, which is posted with the 510(k) boilers on the H drive.]		

For a submission, which relies on a recognized standard without a
declaration of conformity, a statement that the manufacturer
intends to conform to a recognized standard and that supporting
data will be available before marketing the device.
For a submission, which relies on a non-recognized standard that
has been historically accepted by FDA, a statement that the
manufacturer intends to conform to a recognized standard and
that supporting data will be available before marketing the device.
For a submission, which relies on a non-recognized standard that
has not been historically accepted by FDA, a statement that the
manufacturer intends to conform to a recognized standard and
that supporting data will be available before marketing the device
and any additional information requested by the reviewer in order
to determine substantial equivalence.
Any additional information, which is not covered by the guidance
document, special control, recognized standard and/or non-
recognized standard, in order to determine substantial
equivalence.

\* - When completing the review of an abbreviated 510(k), please fill out an Abbreviated Standards Data Form (located on the H drive) and list all the guidance documents, special controls, recognized standards and/or non-recognized standards, which were noted by the sponsor.

Section 4: Additional Requirements for ABBREVIATED and TRADITIONAL 510(k) submissions (If Applicable):

	Present	Inadequate or Missing
a) Biocompatibility data for all patient-contacting materials, OR certification of identical material/formulation:		
b) Sterilization and expiration dating information:	1/	
i) sterilization process  ii) validation method of sterilization process  iii) SAL		
iv) packaoino		
v) specify pyrogen free vi) ETO residues		
vii) radiation dose viii) Traditional Method or Non-Traditional Method c) Software Documentation:		

Items with checks in the "Present or Adequate" column do not require e additional information from the sponsor. Items with checks in the "Missing or Inadequate" column must be submitted before substantive review of the document.

Passed ScreeningYesNo Reviewer:No	
Concurrence by Review Branch:	
Date	

The deficiencies identified above represent the issues that we believe need to be resolved before our review of your 510(k) submission can be successfully completed. In developing the deficiencies, we carefully considered the statutory criteria as defined in Section 513(i) of the Federal Food, Drug, and Cosmetic Act for determining substantial equivalence of your device. We also considered the burden that may be incurred in your attempt to respond to the deficiencies. We believe that we have considered the least burdensome approach to resolving these issues. If, however, you believe that information is being requested that is not relevant to the regulatory decision or that there is a less burdensome way to resolve the issues, you should follow the procedures outlined in the "A Suggested Approach to Resolving Least Burdensome Issues" document. It is available on our Center web page at: http://www.fda.gov/cdrh/modact/leastburdensome.html

THE 510(K) DOCUMENTATION FORMS ARE AVAILABLE ON THE LAN UNDER 510(K) BOILERPLATES TITLED "DOCUMENTATION" AND MUST BE FILLED OUT WITH EVERY FINAL DECISION (SE, NSE, NOT A DEVICE, ETC.).

## "SUBSTANTIAL EQUIVALENCE" (SE) DECISION MAKING DOCUMENTATION

	K	<del></del>	
Revie	wer:		, , , , , , , , , , , , , , , , , , ,
Divis	ion/Branch:	*	
Devic	e Name:		
Produ	ct To Which Compared (510(K) Number If H	Known):	
		YES NO	
1.	Is Product A Device		If NO = Stop
2.	Is Device Subject To 510(k)?		If NO = Stop
3.	Same Indication Statement?		If YES = Go To 5
4.	Do Differences Alter The Effect Or Raise New Issues of Safety Or Effectiveness?		If YES = Stop NE
5.	Same Technological Characteristics?		If YES = Go To 7
6.	Could The New Characteristics Affect Safety Or Effectiveness?		If YES = Go To 8
7.	Descriptive Characteristics Precise Enough?		If NO = Go To 10 If YES = Stop SE
8.	New Types Of Safety Or Effectiveness Questions?		If YES = Stop NE
9.	Accepted Scientific Methods Exist?		If NO = Stop NE
10.	Performance Data Available?		If NO = Request Data
11.	Data Demonstrate Equivalence?		Final Decision:

Note: In addition to completing the form on the LAN, "yes" responses to questions 4, 6, 8, and 11, and every "no" response requires an explanation.

- 1. Intended Use:
- 2. Device Description: Provide a statement of how the device is either similar to and/or different from other marketed devices, plus data (if necessary) to support the statement. Is the device life-supporting or life sustaining? Is the device implanted (short-term or long-term)? Does the device design use software? Is the device sterile? Is the device for single use? Is the device over-the-counter or prescription use? Does the device contain drug or biological product as a component? Is this device a kit? Provide a summary about the devices design, materials, physical properties and toxicology profile if important.

## EXPLANATIONS TO "YES" AND "NO" ANSWERS TO QUESTIONS ON PAGE 1 AS NEEDED

- 1. Explain why not a device:
- Explain why not subject to 510(k):
- 3. How does the new indication differ from the predicate device's indication:
- 4. Explain why there is or is not a new effect or safety or effectiveness issue:
- 5. Describe the new technological characteristics:
- 6. Explain how new characteristics could or could not affect safety or effectiveness:
- 7. Explain how descriptive characteristics are not precise enough:
- 8. Explain new types of safety or effectiveness questions raised or why the questions are not new:
- 9. Explain why existing scientific methods can not be used:
- 10. Explain what performance data is needed:
- 11. Explain how the performance data demonstrates that the device is or is not substantially equivalent:

ATTACH ADDITIONAL SUPPORTING INFORMATION

## 510(k) MEMORANDUM

Date:

12/20/04

Reviewer:

Christopher Hack BSE, Biomedical Engineer

FDA/CDRH/ODE/DGRD/ORDB, HFZ-410

Document #:

K042841

Date on Submission: 10/13/04

Review Initiated:

12/12/04

Received in ODE: 10/14/04

Due Date (60 Days):

Document Received: 10/18/04

Decision Date (75 Days): 12/28/04

**RECOMMENDATION: SE** 

Sponsor and Official Contact:

Patricia S Andborn Beres

56 East Bell Dr. P.O. Box 587 Warsaw, IN 46582

Establishment Registration Number

1825034

## INTERNAL ADMINISTRATIVE FORM

		YES	NO
1.	Did the firm request expedited review?		х
2.	Did we grant expedited review?		X
3.	Have you verified that the Document is labeled Class III for GMP purposes?	Х	
4.	If, not, has POS been notified?		<u> </u>
5.	Is the product a device?	Х	- "
6.	Is the device exempt from 510(k) by regulation or policy?		x
7.	Is the device subject to review by CDRH?	x	
8.	Are you aware that this device has been the subject of a previous NSE decision?		X
9.	If yes, does this new 510(k) address the NSE issue(s), (e.g., performance data)?		
10.	Are you aware of the submitter being the subject of an integrity investigation?		X
	If, yes, consult the ODE Integrity Officer.		
	Has the ODE Integrity Officer given permission to proceed with the review? (Blue Book Memo #I91-2 and Federal Register 90N0332, September 10, 1991.		

## DECISION MAKING FLOWCHART FOR "SUBSTANTIAL EQUIVALENCE"

		YES	NO	
1.	Is the product a device?	X		NO then Stop
2.	Is the device subject to 510(k)?	Х	-	NO then Stop
3.	Is the indication statement the same?	X		YES then Go To 5
4.	Do differences in the indication statement raise new issues of safety and effectiveness?			YES then NSE
5.	Does the device have the same technological characteristics?		X	YES then Go To 7
6.	Could the new characteristics affect safety and effectiveness?			YES then Go To 8
7.	Are the descriptive characteristics precise enough?	X		NO then Go To 10 YES then SE
8.	Are there new types of safety and effectiveness questions?			YES then NSE
9.	Do accepted scientific methods exist to test the impact of the new characteristics?			NO then NSE
	Is performance data available?			NO then Request Data
11.	Does the performance data demonstrate substantial equivalence?			FINAL DECISION: SE

#### **SUMMARY OF REVIEW**

This is a metal on metal acetabular system. Only minor modifications were made in this submission. There is an addition of a slightly rougher surface for the liner cup interface, and there is the addition of two screw holes. Neither of these additions warranted addition testing information, however, the taper roughening was evaluated with push out testing to show that the taper was indeed strengthened by the slight surface roughening. The taper angles remained the same. I recommend that this device be found SE.

#### **REQUIRED FORMS**

Truthful & Accuracy Statement:
Indications for Use Page (IFU):
510k Summary:

#### **DEVICE IDENTIFICATION**

Trade Name:  $M^2a / C^2a$  Acetabular system, metallic acetabular system

Regulation Number: 888.3330 – hip joint metal/semi constrained, with uncemented acetabular component prosthesis

888.3320 – hip joint metal / metal semi-constrained acetabular component

prosthesis 888.3320 – hip joint metal / metal semi-constrained acetabular component

Regulatory Class: Class III Product Code: KWA, JDL

#### INTENDED USE AND INDICATIONS

The M<sup>2</sup>a/C<sup>2</sup>a acetabular system us intended for cemented or non-cemented use in cases of:

- 1) non-inflammatory degenerative joint disease including osteoarthritis and avascular necrosis
- 2) rheumatoid arthritis
- 3) correction of functional deformity
- 4) revision procedures where other treatment or devices have been unsuccessful
- 5) treatment of non-union, femoral neck fracture, trochanteric fractures of the proximal femur with head involvement, unmanageable using other techniques

#### DEVICE DESCRIPTION

The M²a/C²a Acetabular System consists of a titanium outer acetabular shell with a cobalt alloy metallic liner for metal on metal articulation. A cobalt alloy modular head completes the system. The acetabular shells are hemispherical in shape to closely match the natural acetabulum. The shells are available in out diameters of 48mm to 70 mm in 2mm increments. The shell features eight radial fins to aid in the prevention of rotation. Two screw holes in the dome allow for additional fixation by the use of 6.5mm screws. The outer surfaces of the shells are covered with Biomet's plasma spray coating.

The metallica liners contained in this submission are identical to those previously cleared in K993438-metal on metal acetabular system and  $K003363 - M^2a$  32mm taper system. The taper locking mechanism of the shells is also identical to that of the shells cleared in these two 510(k)s.

The metallic cobalt alloy bearing liner fits into the outer shell by means of a taper similar to the taper used for attachment of a modular head to a femoral stem. The locking mechanism consists of an 18°55' taper.

than half the allowable deviation set forth in the ISO standard for modular head sphericity (ISO 7206-2)

The metallic liners articulate with cobalt alloy modular heads identical to those cleared through K993438 an dK0033633- M<sup>2</sup>a 32mm taper system. (b) (4)

Both 28mm and 32mm cobalt alloy modular heads are available in seven neck lengths ranging from -6mm to +12mm. Each modular head has Biomet's type 1 taper and will mate with any Biomet's type 1 taper femoral component.

The acetabular shell diameters and overall geometry are identical to the Mallory/ head shells cleared K993438 and K003363. The only difference is the addition of 2 screw holes to the Mallory/ head configuration. Besides the Mallory head style shells, the predicate 510(k)s contained the universal acetabular shell with 2 screw holes in the dome.

The shells inner taper angle that provides fixation of the metallic liner is identical to the predicates. Slight dimensional changes have been made to the inner shell geometry to insure proper seating of the liner in the shell even when the tolerances are at their limits. The surface roughness of the shell's inner taper has been increased over the predicate.

## **DEVICE MATERIALS**

Acetabular shell is manufactured form ti-alloy (ti-6al-4v) conforming to ASTM F-136 or F-1472

The outer surface of the shell is spray coated with ti-alloy powered conforming to ASTM F-1580. The metallic liner and modular heads are manufactured for cobalt chrome alloy CoCrMo conforming to ASTM F-1537.



## **STERILIZATION**

Radiation type: Gamma

Source: cobalt 60 Min dose; 25Kgy Max dose: 40kGy

Sal 10<sup>-6</sup>

No pyrogen free claims

#### **LABELING**

Adequate

#### **TESTING DETAILS**


Liner locking mechanism: The only modifications the sponsor had made to the liner locking mechanism is by adding a roughened surface finish to the female taper of the Acetabular shell. The sponsor says this modification was done in order to be able to use the same shell with a future ceramic system. The sponsor has been able to show that the push –out strength with the ceramic liners in place was increased 15lbs from 165lbs to 180 lbs with the female taper roughening. The sponsor feels that the same result will be replicated with the metal liners.

Although the sponsor did not testing the roughening with the metal liners, the data from the ceramic liners with the roughened acetabular cup show that the taper is strengthened by the roughening. It does not however show by how much the taper is strengthened.

## ADDITIONAL INFORMATION RELATING TO PREDICATE DEVICE

	New device	Metal on metal Acetabular system		M <sup>2</sup> a 32mm Taper System		Mallory Head PF Acetabular Component
Manufacture	Biomet	Biomet		Biomet		Biomet
K #	042841	K993438		K003363		K861114
Intended use	Cemented and uncemented	Uncemented		Uncemented		Cemented
			Liners			
Liner material	Co-Cr-Mo	Co-Cr-Mo		Co-Cr-Mo		n/a
Liner internal diameter	28mm &32mm	28mm		32mm		n/a
Liner size designation	37mm &41mm	37mm		41mm		n/a
Liner locking mechanism	Taper	Тарег		Taper		n/a
		Ace	etabular Shell			
Shell style	Mallory/head	Mallory/ head	Universal	Mallory/ head	Universal	Mallory / head
Shell material	Ti-6Al-4V	Ti-6Al-4V	Ti-6Al-4V	Ti-6Al-4V	Ti-6AI-4V	Ti-6Al-4V
Liner locking mechanism	Taper	Taper	Taper	Taper	Taper	Hex-loc
Shell outer	48mm-70mm	52mm –	48mm-	52mm-	52mm-	46mm – 70mm
diameters	by 2mm increments	70 mm by 2mm increments	70mm by 2mm increments	70mm by 2mm increments	70mm by 2mm increments	by 4mm increments
Shell features	Radial fins 2 holes plasma spray angled cutouts	Radial fins solid plasma spray straight cut-outs	2 holes plasma spray straight cut outs	Radial fins solid plasma spray straight cut outs	2 holes plasma spray straight cut outs	Fins plasma spray

## 510(k) "SUBSTANTIAL EQUIVALENCE" DECISION-MAKING PROCESS



- 510(k) Submissions compare new devices to marketed devices. FDA requests additional information if the relationship between marketed and "predicate" (pre-Amendments or reclassified post-Amendments) devices is unclear.
- \*\* This decision is normally based on descriptive information alone, but limited testing information is sometimes required.
- ♦ ♦ ♦ Data maybe in the 510(k), other 510(k)s, the Center's classification files, or the literature.