K042911

DEC 27 2004

**Cementek® LV**  
**510(k) Summary**  
**October 18, 2004**

**Submitter:** Teknimed, S.A.  
11 rue Apollo  
31240 L'Union  
FRANCE

**Contact person:** J.D. Webb  
1001 Oakwood Blvd  
Round Rock, TX 78681  
512-388-0199

**Trade Name:** Cementek® LV

**Common name:** Bone void filler

**Classification name:** Class II per 21 CFR section 888.3045

**Product Code:** MQV

**Equivalent Device:** Cementek ® (K040669)

**Device Description**

As an injectable bone substitute, Cementek® LV is packaged as a solid phase and a liquid phase. The liquid and solid phases are mixed in the operating room, then introduced with a syringe into the osseous cavity and allowed to set. This reaction is an athermic reaction resulting in an apatitic calcium phosphate cement. Cementek® LV is marketed in a 16cc dosage.

**Intended Use**

Cementek® LV is intended for use only as a bone void filler for voids or gaps that are not intrinsic to the stability of the bony structure. Cementek® LV is indicated for use in the treatment of surgically created osseous defects or osseous defects created from traumatic injury to bone. The putty can be injected into the bony voids or gaps in the skeletal system (i.e. extremities, spine, and pelvis). Following placement in the bony voids or gap, the calcium phosphate scaffold resorbs and is replaced with bone during the healing process.

**Summary of Technological Characteristics Compared to Predicate Device**

Cementek® LV is equivalent to Cementek® in terms of physical form, how supplied, compressive strength, porosity, average pore size, composition of final product and indications.

DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

DEC 27 2004

Teknimed, S.A.  
c/o Mr. J.D. Webb  
The Orthomedix Group, Inc.  
1001 Oakwood Blvd.  
Round Rock, TX 78681

Re: K042911  
Trade Name: Cementek LV Bone Substitute  
Regulation Number: 21 CFR 888.3045  
Regulation Name: Resorbable calcium salt bone void filler device  
Regulatory Class: Class II  
Product Code: MQV  
Dated: December 6, 2004  
Received: December 8, 2004

Dear Mr. Webb:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good


Page 2 – Mr. J.D. Webb

manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read "Celia M. Witten". The signature is written in a cursive style with a large initial "C" and "W".

Celia M. Witten, Ph.D., M.D.  
Director

Division of General, Restorative and  
Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K042911

Device Name: Cementek® LV

### Indications for Use:

Cementek® LV is intended for use as a bone void filler for voids or gaps that are not intrinsic to the stability of the bony structure. Cementek® LV is indicated for use in the treatment of surgically created osseous defects or osseous defects created from traumatic injury to bone. The putty can be injected into the bony voids or gaps in the skeletal system (i.e. extremities, spine, and pelvis). Following placement in the bony voids or gap, Cementek® LV resorbs and is replaced with bone during the healing process.

Prescription Use  X   
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

*for Mark A. Melhem*  
\_\_\_\_\_  
**(Division Sign-Off)**  
**Division of General, Restorative,**  
**and Neurological Devices**

**510(k) Number** K042911

**DEPARTMENT OF HEALTH & HUMAN SERVICES**

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

DEC 27 2004

Teknimed, S.A.  
c/o Mr. J.D. Webb  
The Orthomedix Group, Inc.  
1001 Oakwood Blvd.  
Round Rock, TX 78681

Re: K042911  
Trade Name: Cementek LV Bone Substitute  
Regulation Number: 21 CFR 888.3045  
Regulation Name: Resorbable calcium salt bone void filler device  
Regulatory Class: Class II  
Product Code: MQV  
Dated: December 6, 2004  
Received: December 8, 2004

Dear Mr. Webb:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good


Page 2 – Mr. J.D. Webb

manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

Celia M. Witten, Ph.D., M.D.  
Director  
Division of General, Restorative and  
Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K042911

Device Name: Cementek® LV

### Indications for Use:

Cementek® LV is intended for use as a bone void filler for voids or gaps that are not intrinsic to the stability of the bony structure. Cementek® LV is indicated for use in the treatment of surgically created osseous defects or osseous defects created from traumatic injury to bone. The putty can be injected into the bony voids or gaps in the skeletal system (i.e. extremities, spine, and pelvis). Following placement in the bony voids or gap, Cementek® LV resorbs and is replaced with bone during the healing process.

Prescription Use X  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE  
OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

*for Mark A. Millheim*  
**(Division Sign-Off)**  
**Division of General, Restorative,  
and Neurological Devices**

510(k) Number K042911

**DEPARTMENT OF HEALTH & HUMAN SERVICES**

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

**NOV 15 2004**

Teknimed, S.A.  
c/o Mr. J.D. Webb  
The Orthomedix Group, Inc.  
1001 Oakwood Blvd.  
Round Rock, TX 78681

Re: K042911  
Trade Name: Cementek LV Bone Substitute  
Dated: October 18, 2004  
Received: October 21, 2004

Dear Mr. Webb:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above. We cannot determine if the device is substantially equivalent to a legally marketed predicate device based solely on the information you provided. To complete the review of your submission, we require the following additional information.

(b)(4)

Handwritten initials or signature in the bottom right corner of the page.

Page 2 - Mr. J.D. Webb

(b)(4)

The deficiencies identified above represent the issues that we believe need to be resolved before our review of your 510(k) submission can be successfully completed. In developing the deficiencies, we carefully considered the statutory criteria as defined in Section 513(i) of the Federal Food, Drug, and Cosmetic Act for determining substantial equivalence of your device.

We also considered the burden that may be incurred in your attempt to respond to the deficiencies. We believe that we have considered the least burdensome approach to resolving these issues. If, however, you believe that information is being requested that is not relevant to the regulatory decision or that there is a less burdensome way to resolve the issues, you should follow the procedures outlined in the "A Suggested Approach to Resolving Least Burdensome Issues" document. It is available on our Center web page at: <http://www.fda.gov/cdrh/modact/leastburdensome.html>

You may not market this device until you have provided adequate information described above and required by 21 CFR 807.87(l), and you have received a letter from FDA allowing you to do so. If you market the device without conforming to these requirements, you will be in violation of the Federal Food, Drug, and Cosmetic Act (Act). You may, however, distribute this device for investigational purposes to obtain clinical data if needed to establish substantial equivalence. Clinical investigations of this device must be conducted in accordance with the investigational device exemption (IDE) regulations.

If the information, or a request for an extension of time, is not received within 30 days, we will consider your premarket notification to be withdrawn and your submission will be deleted from our system. If you submit the requested information after 30 days it will be considered and processed as a new 510(k); therefore, all information previously submitted must be resubmitted so that your new 510(k) is complete. Please note our guidance document entitled, "Guidance for Industry and FDA Staff FDA and Industry Actions on Premarket Notification (510(k)) Submissions: Effect on FDA Review Clock and Performance Assessment". The purpose of this document is to assist agency staff and the device industry in understanding how various FDA and industry actions that may be taken on 510(k)s should affect the review clock for purposes of meeting the Medical Device User Fee and Modernization Act. You may review this document at <http://www.fda.gov/cdrh/mdufma/guidance/1219.html>.

18


Page 3 - Mr. J.D. Webb

The requested information, or a request for an extension of time, should reference your above 510(k) number and should be submitted in duplicate to:

Food and Drug Administration  
Center for Devices and  
Radiological Health  
Document Mail Center (HFZ-401)  
9200 Corporate Boulevard  
Rockville, Maryland 20850

If you have any questions concerning the contents of the letter, please contact Nadine Y. Sloan at (301)594-1296, ext 205. If you need information or assistance concerning the IDE regulations, please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,

For Celia M. Witten, Ph.D., M.D.

Director  
Division of General, Restorative and  
Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

Teknimed, S.A.  
c/o Mr. J.D. Webb  
The Orthomedix Group, Inc.  
1001 Oakwood Blvd.  
Round Rock, TX 78681

Re: K042911  
Trade Name: Cementek LV Bone Substitute  
Dated: October 18, 2004  
Received: October 21, 2004

Dear Mr. Webb:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above. We cannot determine if the device is substantially equivalent to a legally marketed predicate device based solely on the information you provided. To complete the review of your submission, we require the following additional information.

(b)(4)

Page 2 - Mr. J.D. Webb

(b)(4)

The deficiencies identified above represent the issues that we believe need to be resolved before our review of your 510(k) submission can be successfully completed. In developing the deficiencies, we carefully considered the statutory criteria as defined in Section 513(i) of the Federal Food, Drug, and Cosmetic Act for determining substantial equivalence of your device.

We also considered the burden that may be incurred in your attempt to respond to the deficiencies. We believe that we have considered the least burdensome approach to resolving these issues. If, however, you believe that information is being requested that is not relevant to the regulatory decision or that there is a less burdensome way to resolve the issues, you should follow the procedures outlined in the "A Suggested Approach to Resolving Least Burdensome Issues" document. It is available on our Center web page at: <http://www.fda.gov/cdrh/modact/leastburdensome.html>

You may not market this device until you have provided adequate information described above and required by 21 CFR 807.87(l), and you have received a letter from FDA allowing you to do so. If you market the device without conforming to these requirements, you will be in violation of the Federal Food, Drug, and Cosmetic Act (Act). You may, however, distribute this device for investigational purposes to obtain clinical data if needed to establish substantial equivalence. Clinical investigations of this device must be conducted in accordance with the investigational device exemption (IDE) regulations.

If the information, or a request for an extension of time, is not received within 30 days, we will consider your premarket notification to be withdrawn and your submission will be deleted from our system. If you submit the requested information after 30 days it will be considered and processed as a new 510(k); therefore, all information previously submitted must be resubmitted so that your new 510(k) is complete. Please note our guidance document entitled, "Guidance for Industry and FDA Staff FDA and Industry Actions on Premarket Notification (510(k)) Submissions: Effect on FDA Review Clock and Performance Assessment". The purpose of this document is to assist agency staff and the device industry in understanding how various FDA and industry actions that may be taken on 510(k)s should affect the review clock for purposes of meeting the Medical Device User Fee and Modernization Act. You may review this document at <http://www.fda.gov/cdrh/mdufma/guidance/1219.html>.

Page 3 - Mr. J.D. Webb

The requested information, or a request for an extension of time, should reference your above 510(k) number and should be submitted in duplicate to:

Food and Drug Administration  
 Center for Devices and  
 Radiological Health  
 Document Mail Center (HFZ-401)  
 9200 Corporate Boulevard  
 Rockville, Maryland 20850

If you have any questions concerning the contents of the letter, please contact Nadine Y. Sloan at (301)594-1296, ext 205. If you need information or assistance concerning the IDE regulations, please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,

Celia M. Witten, Ph.D., M.D.  
 Director  
 Division of General, Restorative and  
 Neurological Devices  
 Office of Device Evaluation  
 Center for Devices and  
 Radiological Health

Enclosure

22

FILE  
 COPY

OFFICE	SURNAME	DATE	OFFICE	SURNAME	DATE	OFFICE	SURNAME	DATE
2410	Sloan	11/15/04						
2410	Stevens	11/15/04						

cc: HFZ-401 DMC  
HFZ-404 510(k) Staff  
HFZ- 410 DGRND  
D.O.

nyr:11-15-04

DEPARTMENT OF HEALTH AND HUMAN SERVICES

Public Health Service

Food and Drug Administration  
Center for Devices and  
Radiological Health  
Office of Device Evaluation  
Document Mail Center (HFZ-401)  
9200 Corporate Blvd.  
Rockville, Maryland 20850

October 21, 2004

TEKNIMED SA  
C/O THE ORTHOMEDIX GROUP, INC.  
1001 OAKWOOD BLVD  
ROUND ROCK, TX 78681  
ATTN: J.D. WEBB

510(k) Number: K042911  
Received: 21-OCT-2004  
Product: CEMENTEK LV

The Food and Drug Administration (FDA), Center for Devices and Radiological Health (CDRH), has received the Premarket Notification you submitted in accordance with Section 510(k) of the Federal Food, Drug, and Cosmetic Act (Act) for the above referenced product. We have assigned your submission a unique 510(k) number that is cited above. Please refer prominently to this 510(k) number in any future correspondence that relates to this submission. We will notify you when the processing of your premarket notification has been completed or if any additional information is required. YOU MAY NOT PLACE THIS DEVICE INTO COMMERCIAL DISTRIBUTION UNTIL YOU RECEIVE A LETTER FROM FDA ALLOWING YOU TO DO SO.

On May 21, 2004, FDA issued a Guidance for Industry and FDA Staff entitled, "FDA and Industry Actions on Premarket Notification (510(k)) Submissions: Effect on FDA Review Clock and Performance Assessment". The purpose of this document is to assist agency staff and the device industry in understanding how various FDA and industry actions that may be taken on 510(k)s should affect the review clock for purposes of meeting the Medical Device User Fee and Modernization Act. Please review this document at <http://www.fda.gov/cdrh/mdufma/guidance/1219.html>.

Please remember that all correspondence concerning your submission MUST be sent to the Document Mail Center (DMC)(HFZ-401) at the above letterhead address. Correspondence sent to any address other than the one above will not be considered as part of your official premarket notification submission. Also, please note the new Blue Book Memorandum regarding Fax and E-mail Policy entitled, "Fax and E-Mail Communication with Industry about Premarket Files Under Review". Please refer to this guidance for information on current fax and e-mail practices at [www.fda.gov/cdrh/ode/a02-01.html](http://www.fda.gov/cdrh/ode/a02-01.html).

You should be familiar with the regulatory requirements for medical device available at Device Advice <http://www.fda.gov/cdrh/devadvice/>. If you have other procedural or policy questions, or want information on how to check on the status of your submission, please contact DSMICA at (301) 443-6597 or its toll-free number (800) 638-2041, or at their Internet address <http://www.fda.gov/cdrh/dsmamain.html> or me at (301)594-1190.

Sincerely yours,

Marjorie Shulman  
Supervisory Consumer Safety Officer  
Office of Device Evaluation  
Center for Devices and Radiological Health

K042911

DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION  <b>MEDICAL DEVICE USER FEE COVER SHEET</b>	PAYMENT IDENTIFICATION NUMBER (b) (4) Write the Payment Identification Number on your check.
---	---

A completed Cover Sheet must accompany each original application or supplement subject to fees. The following actions must be taken to properly submit your application and fee payment:

1. Electronically submit the completed Cover Sheet to the Food and Drug Administration (FDA) before payment is sent.
2. Include a printed copy of this completed Cover Sheet with a check made payable to the Food and Drug Administration. Remember that the Payment Identification Number must be written on the check.
3. Mail Check and Cover Sheet to the US Bank Lock Box, FDA Account, P.O. Box 956733, St. Louis, MO 63195-6733. (Note: In no case should payment be submitted with the application.)
4. If you prefer to send a check by a courier, the courier may deliver the check and Cover Sheet to: US Bank, Attn: Government Lockbox 956733, 1005 Convention Plaza, St. Louis, MO 63101. (Note: This address is for courier delivery only. Contact the US Bank at 314-418-4821 if you have any questions concerning courier delivery.)
5. For Wire Transfer Payment Procedures, please refer to the MDUFMA Fee Payment Instructions at the following URL: <http://www.fda.gov/cdrh/mdufma/faqs.html#3a>. You are responsible for paying all fees associated with wire transfers.
6. Include a copy of the completed Cover Sheet in volume one of the application when submitting to the FDA at either the CBER or CDRH Document Mail Center.

1. COMPANY NAME AND ADDRESS (Include name, street address, city, state, country, and post office code)  TEKNIMED 11 RUE APOLLO L1 UNION, 31240 FR  1.1 EMPLOYER IDENTIFICATION NUMBER (EIN)	2. CONTACT NAME JD WEBB  2.1 E-MAIL ADDRESS ortho.medix@sbcglobal.net  2.2 TELEPHONE NUMBER (Include Area Code) 512-388-0199  2.3 FACSIMILE (FAX) NUMBER (Include Area Code) 512-388-0199
--	---

3. TYPE OF PREMARKET APPLICATION (Select one of the following in each column; if you are unsure, please refer to the application descriptions at the following web site: <http://www.fda.gov/oc/mdufma>)

Select an application type: Premarket notification (510(k)); except for third party Biologic License Application (BLA) Premarket Approval Application (PMA) Modular PMA Product Development Protocol (PDP) Premarket Report (PMR)	3.1 Select one of the types below: Original Application  Supplement Types: Efficacy (BLA, PMR) Panel Track (PMA, PMR, PDP) Real-Time (PMA, PMR, PDP) 180-day (PMA, PMR, PDP)
---	---

4. ARE YOU A SMALL BUSINESS? (See the instructions for more information on determining this status.)

YES, I meet the small business criteria and have submitted the required qualifying documents to FDA      NO, I am not a small business

4.1 If Yes, please enter your Small Business Decision Number:

5. IS THIS PREMARKET APPLICATION COVERED BY ANY OF THE FOLLOWING USER FEE EXCEPTIONS? IF SO, CHECK THE APPLICABLE EXCEPTION.

This application is the first PMA submitted by a qualified small business, including any affiliates, parents, and partner firms	The sole purpose of the application is to support conditions of use for a pediatric population
This biologic application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only	The application is submitted by a state or federal government entity for a device that is not to be distributed commercially

6. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION FOR WHICH FEES WERE WAIVED DUE TO SOLE USE IN A PEDIATRIC POPULATION THAT NOW PROPOSES CONDITION OF THE USE FOR ANY ADULT POPULATION? (If so, the application is subject to the fee that applies for an original premarket approval application (PMA).)

YES      NO

7. USER FEE PAYMENT AMOUNT SUBMITTED FOR THIS APPLICATION

(b) (4)

SK32  
 PM  
 II 30

# THE ORTHOMEDIX GROUP, INC.

October 18, 2004

Document Mail Center (HFZ401)  
Center for Devices and Radiological Health  
Food and Drug Administration  
9200 Corporate Boulevard  
Rockville, Maryland 20850

RECEIVED  
OCT 21 2004  
CDRH

Re: Premarket Notification – Cemetek® LV Bone Substitute  
Modification to Cemetek® Bone Substitute K041493


Dear Sir/Madam,

Teknimed, S.A. hereby submits this **Special 510(k): Device Modification** to request modifications to our Cemetek® Bone Substitute (K041493). We believe the modifications to this product are eligible for the Special 510(k) process since they have the same fundamental scientific technology and intended use as the predicate device.

Applicant	Teknimed, S.A.
Address	11 rue Apollo 31240 L'Union FRANCE
Contact Person	J.D. Webb
Address	1001 Oakwood Blvd Round Rock, TX 78681
Telephone Number	512-388-0199
Fax Number	512-388-0199
E-mail Address	ortho.medix@sbcglobal.net

We consider our intent to market this device as confidential commercial information and request that it be treated as such by FDA. We have taken precautions to protect the confidentiality of the intent to market this device.

Respectfully,

J.D. Webb  
Authorized Contact Person

SK32


## **Table of Contents**

Truth and Accurate Statement .....	Page ii
Indications for Use.....	Page iii
510(k) Summary .....	Page iv
Device Name.....	Page 1
Registration Number.....	Page 1
Classification, Panel & Product Code.....	Page 1
Performance Standards .....	Page 1
Labeling .....	Page 1
Substantial Equivalence Comparison.....	Page 1
Device Description .....	Page 1
Design Control Activities .....	Page 5
Sterility.....	Page 6
Sample Labeling .....	Exhibit I
Equivalent Device.....	Exhibit II
Manufacturing Process .....	Exhibit III
Test Report.....	Exhibit IV
Risk Analysis .....	Exhibit V

**PREMARKET NOTIFICATION**  
**TRUTHFUL AND ACCURATE STATEMENT**

[As required by 21 CFR 807.87(k)]

I certify that, in my capacity as *<M.LEONARD Alain >* of *<TEKNIMED>* I believe to the best of my knowledge, that all data and information submitted in the premarket notification are truthful and accurate and that no material fact has been omitted.

  
\_\_\_\_\_  
(Signature)

*A LEONARD*  
\_\_\_\_\_  
(Typed name)

*06.10.2004*  
\_\_\_\_\_  
(Dated)

\_\_\_\_\_  
\*(Premarket Notification [510(k)] Number)

\* For a new submission, leave the space for the 510(k) number blank.

## Indications for Use

510(k) Number (if known): K042911

Device Name: Cementek® LV

### Indications for Use:

Cementek® LV is intended for use as a bone void filler for voids or gaps that are not intrinsic to the stability of the bony structure. Cementek® LV is indicated for use in the treatment of surgically created osseous defects or osseous defects created from traumatic injury to bone. The putty can be injected into the bony voids or gaps in the skeletal system (i.e. extremities, spine, and pelvis). Following placement in the bony voids or gap, Cementek® LV resorbs and is replaced with bone during the healing process.

Prescription Use  X  AND/OR Over-The-Counter Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D) (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE  
OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

**Cementek® LV**  
**510(k) Summary**  
**October 18, 2004**

**Submitter:** Teknimed, S.A.  
11 rue Apollo  
31240 L'Union  
FRANCE

**Contact person:** J.D. Webb  
1001 Oakwood Blvd  
Round Rock, TX 78681  
512-388-0199

**Trade Name:** Cementek® LV

**Common name:** Bone void filler

**Classification name:** Class II per 21 CFR section 888.3045

**Product Code:** MQV

**Equivalent Device:** Cementek ® (K040669)

**Device Description**

As an injectable bone substitute, Cementek® LV is packaged as a solid phase and a liquid phase. The liquid and solid phases are mixed in the operating room, then introduced with a syringe into the osseous cavity and allowed to set. This reaction is an athermic reaction resulting in an apatitic calcium phosphate cement. Cementek® LV is marketed in a 16cc dosage.

**Intended Use**

Cementek® LV is intended for use only as a bone void filler for voids or gaps that are not intrinsic to the stability of the bony structure. Cementek® LV is indicated for use in the treatment of surgically created osseous defects or osseous defects created from traumatic injury to bone. The putty can be injected into the bony voids or gaps in the skeletal system (i.e. extremities, spine, and pelvis). Following placement in the bony voids or gap, the calcium phosphate scaffold resorbs and is replaced with bone during the healing process.

**Summary of Technological Characteristics Compared to Predicate Device**

Cementek® LV is equivalent to Cementek® in terms of physical form, how supplied, compressive strength, porosity, average pore size, composition of final product and indications.

October 14, 2004  
Teknimed S.A.  
Cementek® LV  
Special 510(k)  
Page 1 of 6

---

**1. Device Name**

- 1.1. Trade Name  
Cementek® LV
- 1.2. Common Name  
Bone void filler
- 1.3. Classification Name  
Filler, calcium sulfate, preformed pellets

**2. Registration Number**

9615788

**3. Classification**

- 3.1. Class  
Class II per 21 CFR Sec. 888.3045
- 3.2. Panel  
Orthopedic
- 3.3. Product Code  
MQV

**4. Performance Standards**

There are no performance standards for calcium phosphate cement.

**5. Labeling (Exhibit I)**

Draft labels and Instructions for Use can be found in Exhibit I.  
No changes to the labels or Instructions for Use have occurred.

**6. Substantial Equivalence Comparison**

Cementek® LV is similar to Cementek® (K041493). A table comparing the various properties is shown in Exhibit II.

**7. Device Description**

As an injectable bone substitute, Cementek® LV is packaged as a solid phase and a liquid phase. The liquid and solid phases are mixed in the operating room, then injected as a paste into the osseous cavity and allowed to set.

7.1. Indications for Use

Cementek® LV is intended for use only as a bone void filler for voids or gaps that are not intrinsic to the stability of the bony structure. Cementek® LV is indicated

October 14, 2004  
Teknimed S.A.  
Cementek® LV  
Special 510(k)  
Page 2 of 6

---

for use in the treatment of surgically created osseous defects or osseous defects created from traumatic injury to bone. The putty can be injected into the bony voids or gaps in the skeletal system (i.e. extremities, spine, and pelvis). Following placement in the bony voids or gap, the calcium phosphate scaffold resorbs and is replaced with bone during the healing process.


#### 7.2. Chemical Composition

Cementek® LV is manufactured by processing and mixing calcium carbonate and calcium dihydrogenophosphate, creating tetracalcium phosphate(TTCP). In parallel, amorphous tricalcium phosphate(TCP) is processed creating  $\alpha$ -tricalcium phosphate. These two components are mixed with sodium glycerophosphate and polydimethyl siloxane with a solvent (cyclohexane) to produce the solid phase of the product. Following mixing the solid phase is dried to remove the solvent.

A schematic of the manufacturing process can be found in Exhibit III.

Calcium hydroxide, orthophosphoric acid and water are mixed and processed to create the liquid phase of Cementek® LV.

The solid and liquid phases are packaged separately in a pouch and flask then boxed together as a kit. A mixing kit comprised of a spatula, a funnel, a mixer and a plunger is provided. (see photo below) The mixed material is introduced to the surgical site with a syringe.

Spatula, funnel, mixer and plunger  
for mixing Cementek® LV

October 14, 2004  
Teknimed S.A.  
Cemetek<sup>®</sup> LV  
Special 510(k)  
Page 3 of 6

---

7.2.1. Raw materials

(b) (4)

7

7.3.1. Physical form

Injectable paste intended to set *in-vivo*.

7.3.2. Cemetek<sup>®</sup> LV is marketed in 16cc dosage.

(b) (4)

October 14, 2004  
Teknimed S.A.  
Cemetek® LV  
Special 510(k)  
Page 4 of 6

---

(b) (4)

[Redacted]

(b) (4)

[Redacted]

[Redacted]

[Redacted]

[Redacted]

[Redacted]

[Redacted]

October 14, 2004  
Teknimed S.A.  
Cemetek® LV  
Special 510(k)  
Page 5 of 6

---

(b) (4)

A large rectangular area of the document is redacted with a solid grey fill. The text "(b) (4)" is positioned at the top left corner of this redacted area.

#### 7.4.6. Animal study

(b) (4)

A large rectangular area of the document is redacted with a solid grey fill. The text "(b) (4)" is positioned at the top left corner of this redacted area.

#### 7.4.7. Cadaveric study

(b) (4)

A large rectangular area of the document is redacted with a solid grey fill. The text "(b) (4)" is positioned at the top left corner of this redacted area.

### 8. Summary of Design Control Activities

The risk analysis method used to assess the impact of the modifications was a Failure Modes and Effects Analysis (FMEA). Differences between Cemetek® LV and Cemetek®, and the risk analysis are presented in Exhibit V.

October 14, 2004  
Teknimed S.A.  
Cemetek® LV  
Special 510(k)  
Page 6 of 6

---

## 9. Sterility

(b) (4)

### 9.5. Packaging

#### Cementek®LV

The liquid is contained in a 10 ml glass flask.

The solid phase is contained in an aluminum heat sealed pouch.

The both products are packaged in a blister.

This blister is packaged in a tyvek pouch.

#### Mixing tools

Double packaged in tyvek pouches

### 9.6. Pyrogens

This material is not tested for the presence of pyrogens.

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**Powder LV**      2004-07

**Powder for ionic cement mix with liquid**      2009-08

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Tél. 33 (0)5 62 96 88 38

Fax 33 (0)5 62 96 28 72

**Caution: U.S. Federal law restricts this device to sale by or on the order of the physician (or properly licensed practitioner)**

Before using Teknimed products, the operating surgeon should study thoroughly the safety information specified in these instructions as well as the product-specific information (product description, surgical procedures, brochures etc.). The related documentation is obtainable from the respective national representative. At the same time the operating surgeon must be aware of the residual risks associated with the use of the intended products.

**GENERAL INSTRUCTIONS**

The implantation of Teknimed products may only be carried out by surgical staff who possess a thorough knowledge of and experience in the area of joint replacement and, in particular, have mastered the product specific surgical techniques relating to Teknimed products. The particular surgical techniques required for Teknimed bone substitute can be learned at Teknimed distributors.

Cementek® LV bone substitute is in the form of an apatitic paste designed for the osteoconductive replacement of bone.

**COMPOSITION**

This product of synthesis is controlled during all its manufacture, from raw material to the final product.

Composition:

Powder:

tetracalcium phosphate	49%
α tricalcium phosphate	38%
Sodium glycerophosphate	12%
Polydimethylsiloxane	1%

Liquid:

calcium hydroxide	3,4%
phosphoric acid	13,8%
water	82,8%

**PREPARATION FOR USE**

Mix in a cup the whole of the powder and liquid to achieve a homogeneous mixture. The paste obtained will become less and less malleable, until becoming hard after 10 min. Full hardening follows in-situ within 48 to 72 hours. Match the quantity of Cementek® LV to the site of bone defect to obtain the fullest contact with the lost bone.

**INSTRUCTIONS FOR USE**

1. Open the sachet and pour the powder in the mixer.
2. Pour the liquid on the powder.
3. Mix energetically with the mixer the powder and the liquid during 1 to 2 min.
4. Place mixed paste into syringe.
5. After having dried the cavity, put the paste in the bone defect. It is possible to clean and to dry the cavity with STERILE hydrogen peroxide, unless otherwise exceptions.

**PRECAUTIONS FOR USE**

The cavity must be carefully irrigated and dried before application of Cementek® LV.

Cementek® LV has stability of shape after 10 minutes. It is recommended to limit the load applied during the first 48 hours, on large bone defects. Cementek® LV is not intended to provide load-bearing structural support during the healing process, therefore, Cementek® LV is contraindicated where the device is intended as load-bearing structural support in the skeletal system.

Manipulations of the paste should be avoided as it may change the mechanical properties.

**INDICATIONS**

Cementek® LV is intended for use as a bone void filler for voids or gaps that are not intrinsic to the stability of the bony structure. Cementek® LV is indicated for use in the treatment of surgically created osseous defects or osseous defects created from traumatic injury to bone. The putty can be injected into the bony voids or gaps in the skeletal system (i.e. extremities, spine, and pelvis). Following placement in the bony voids or gap, Cementek® LV resorbs and is replaced with bone during the healing process

**CONTRAINDICATIONS**

The same as applying to all bone grafts

- metabolic conditions.
- use in an infected area (osteomyelitis, tuberculosis).
- in an area having no possibility of regeneration or infected bone (risk of sequestration).

**POSSIBLE ADVERSE EFFECTS**

Possible adverse effects include, but are not limited to:

- Wound complications including hematoma, site drainage, bone fracture, infection, and other complications that are possible with any surgery
- Fracture or extrusion of Cementek® LV with or without generation of particulate debris
- Deformity of the bone at the site
- Incomplete, or lack of, osseous ingrowth into bone void, as is possible with any bone void filler.

\*

**STERILIZATION**

Cementek® LV is sterilized with 25 to 40 kGy (2,5 to 4,0 Mrad) gamma radiation. All sterile implants are to be kept unopened in the original packaging until the time of implantation. Before using the implant, the protective packaging should be checked for damage as this could be detrimental to the sterility. Aseptic procedures should be observed when removing the implant from the protective package.

Cementek® LV is delivered sterile under double wrapping, ready for use in an operating room.

**Remarks:**

Verify the integrity of the packaging before use, the guarantee of sterilization is 5 years from the date of sterilization.

Use after the peremptory date is not allowed. All re-sterilization of the product is forbidden use only once.

**PATIENT INFORMATION**

The doctor must draw the patient's attention to the contents of the indications and contraindications paragraph, as well as factors which can impair the results of the operation and to possible complications which can arise as a result of an indication. The patient must also be informed about the measures which the doctor will use to minimize the possible effects of these factors.

**PACKAGING AND STORAGE**

Individually wrapped in quantities.

<u>Designation</u>	<u>Reference</u>	<u>Volume cm<sup>3</sup></u>
Cementek® LV	T8150LV	16cm <sup>3</sup>


Store at room temperature

**Equivalence of Cementek® LV with the Predicate Cementek®**  
Changes are noted in **bold and underlined**

<b>Cementek® LV</b>	<i>Property</i>	<b>Cementek®</b>
---------------------	-----------------	------------------

**Chemical Composition**

(b) (4)

<b>Performance Testing - Bench</b>		
------------------------------------	--	--

(b) (4)

<b>Cementek® LV</b>	<b><i>Property</i></b>	<b>Cementek®</b>
---------------------	------------------------	------------------

(b) (4)

*A*

(b) (4)

## **PROCESS OF CEMENTEK<sup>®</sup> LV**

(b) (4)

(b) (4)

Key Engineering Materials Vols. 192-195 (2001) pp. 789-792

Proceedings of the 13th Int. Symp. on Ceramics in Medicine, Bologna, Italy, 22-26 Nov. (2000) pp. 789-792

## Formulation of an Injectable Phosphocalcium Cement

S. Gonçalves<sup>1,3</sup>, A. Brouchet<sup>4</sup>, M. Frèche<sup>1</sup>, F. Rodriguez<sup>2</sup>, B. Delisle<sup>4</sup>  
and J.L. Lacout<sup>1</sup>

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**Keywords:** Bone Ingrowth, Calcium Phosphate Cement, Polyglycol, Silicon

**Abstract:** In orthopedic surgery, the loss or the reinforcement of osseous substance often requires filling of the defective part. In order to make the surgical operations easier we sought to make an injectable form. This study examined the effect of silicone and polyglycol on the injectability, setting time and mechanical properties of the cement.

The basic solid phase was composed of a mixture of tetracalcium phosphate ( $\text{Ca}_4(\text{PO}_4)_2\text{O}$ ),  $\alpha$ -tricalcium phosphate ( $\text{Ca}_3(\text{PO}_4)_2$ ) and sodium glycerophosphate. The basic liquid phase was made up of lime, orthophosphoric acid and water. Silicone was previously dissolved in cyclohexane and introduced in the solid phase. Polyglycol is a water-soluble compound so it is introduced in the liquid phase.

For the mechanical properties, the strong increase in the percentage of additives decreased the compressive strength. Silicone and polyglycol made it possible to improve viscosity without modifying the basic setting time. The rate of evolution was different with the two different additives. From the data it was possible to optimize the formulation of cements to give predicted properties. Testing the *in vivo* implantation of the cement has already started. Preliminary results show the perfect osteointegration of the new cements without reactions to the foreign body in spite of the presence of silicone.

### Introduction:

Nowadays calcium phosphate cements are being increasingly used in surgery. Their current applications are multiple. However new applications using techniques, such as catheterisation [1,2], are possible but require the use of injectable cement. The osseous catheterisation makes it possible to treat the lesions indirectly decreasing the post operational traumatism and the duration of the periods of convalescence. The aim of this work was to modify the formulation of an already existing calcium phosphate cement Cementek<sup>®</sup>[3] to modify its viscosity while preserving its main properties. Cementek<sup>®</sup> is mainly used for osseous fillings: cotyle filling, stabilisation of the implants without PMMA cement and hip arthroplasties. The global chemical reaction is the following :

with  $a + 3b + 4c = 10$  and  $d + 2b + 2c = 6$ .

The cement results from acido-basic reactions and a dissolution-precipitation process. Its evolution towards apatite, successively passing in the brushite then octocalcium phosphate phases, is now well-known [4,5]. The injection of cement carried out just after the mixing will thus be mainly made up of brushite. The role of the additive will be to improve the rheology of this phase of the cement. The effect of various additives was studied using an experimental plan: silicone and polyglycol seemed to be the best.

### Materials and methods:

The basic cement is composed of two phases:

- a liquid phase containing lime, phosphoric acid and water for injectable preparation.
- a solid phase with tetracalcium and tricalcium phosphates and sodium glycerophosphate.

The two additives are not added cement in the same way. Silicone is dissolved beforehand into cyclohexane then introduced into the solid phase. The solvent is completely evaporated

after 48 hours at ambient temperature. The polyglycol is a water-soluble compound and is therefore introduced in the liquid phase. The cement is presented in a package containing the exact quantities of solid and liquid phases for one dose. The package is sterilized by gamma radiation.

All the different cements were tested for injectability, setting time and mechanical properties. Other tests of standard characterization included infrared spectra (Perkin-Elmer FTIR), X-ray diffraction (CPS 120 INEL), scanning electronic microscopy (SEM: LEO 435 VPS). The compression resistance tests were carried out with samples of standard size (diameter and height), using a Hounsfield Serie S device. The setting time, which corresponds to a penetration resistance of 300 g/mm<sup>2</sup>, was assessed with a penetrometer (TANT2). Injectability was measured using a system developed in the laboratory [6], making it possible to determine an injection pressure. It includes a manometer, a syringe, a catheter and a Tee connector.

The study was undertaken on 19 White New Zealand rabbits (male weighing 4.5 to 5.2 kg), divided into 2 groups of 5 individuals [7,8]. Each group corresponded to a time of given retreat (4 and 11 weeks). Each rabbit remained in an individual cage with free access to water and food. Each received 2 implants: in the femur, cement of reference (A) and injectable cement (B). Both femurs were operated in same operational time by external surgical access of the side condyle under general anaesthesia. A window of 6 mm by 6 was cut from the bone and separated from the spongy subjacent tissue, followed creation of an osseous defect 4.5 mm in diameter and 6 mm length. The osseous cavity was drained and filled manually by cement then closed again by the cortical shutter.

#### Results and discussion

Two experimental plans were used to optimize the percentage of additive. Thus the ideal percentage for silicone is about 1% (table 1). In the other plan, it is necessary to add 2% polyglycol (table 2). Thereafter, the influence of each additive was studied in larger quantities. Thus the percentage of silicone was varied up to 10% and the percentage of polyglycol up to 5%.

		Levels					Step
Coded variable	X1, X2, X3	-1.682	-1	0	1	1.682	$\Delta X_i$
Natural variable	X1= % EL	0	0.406	1	1.594	2	0.594
xi	X2= %HS 15	1	1.406	2	2.594	3	0.594
	X3= %NaGP	8	9.014	10.5	11.986	13	1.486

Table 1: Working range and variation step

		Levels					Step
Coded variable	X1, X2, X3	-1.682	-1	0	1	1.682	$\Delta X_i$
Natural variable	X1= %Si V50	0.5	0.906	1.5	2.094	2.5	0.594
xi	X2= %NaGP	7	8.22	10	11.78	13	1.78
	X3= %H <sub>3</sub> PO <sub>4</sub>	11.5	11.81	12.25	12.69	13	0.44

Table 2: Working range and variation step

Influence of percentage of silicone and polyglycol on the mechanical properties: The influence of silicone on the compressive strength (table 3) is very limited up to 5%. For higher values, the mechanical properties are seriously degraded: in large amounts, the silicone

decreases the mechanical resistance due to the increase of the porosity of cement. In the cement with 10 % silicone, the porosity increased by 47 % to 55 %. For under 3 % polyglycol, the mechanical properties of the cement (table 4) were not modified because the porosity remained unchanged, equal to 52 %. Over 5 % additive, the compressive strength strongly decreased. So polyglycol and silicone can present an interest for quantities added lower than or equal to 3 %.

% SIL	0	0.5	1	3	5	10
Compressive strength (MPa)	12.1	11.7	11.7	11.7	11.8	8.7
Setting time (min)	26	24	24	26	28	29

Table 3: Compressive strength & setting time for cement with silicone

% PEG	0	1	2	3	5
Compressive strength (MPa)	12.1	14.3	11.6	12.0	7.9
Setting time (min)	26	24	24	25	26

Table 4: Compressive strength & setting time for cement with polyglycol

Fig 1: X-ray diffraction diagrams of cement with 0% to 5% of polyglycol


Fig 3: After 4 weeks, periphery of cement colonized by new osseous spans

Fig 4: New osseous spans broadsides by a base of osteoblastic cells (high magnification of fig 3)


Fig 2: X-ray diffraction diagrams of cement with 0% to 10 % of silicone

Influence of percentage of silicone and polyglycol on evolution: The figures present the evolution of cements according to the percentages of added additives. The X-ray diffraction diagrams with polyglycol (fig 1) are identical, so the influence on the evolution of cement is limited. However the diagrams (fig 2) with silicone, are different. The rate of evolution of the hydroxyapatite is slower in this case. The major  $\beta$ -TCP phase on the diagrams with 10 % silicone confirms this.

Influence of percentage of silicone and polyglycol on setting time: All the setting times are close to the value of the initial cement (table 3 and 4). The influence of silicone and polyglycol on setting time is limited.

Influence of percentage of silicone and polyglycol on injectability: All cements resulting from the 2 experimental plans were injectable over about thirty centimetres through catheters of 2 or 3 mm internal diameter. These cements have pressures of injection lower than 1 bar and do not surbed in water. The time of injection is related to the percentage of additives but also to the phosphoric acid percentage. All times of injection were lower than 15 minutes.

Animal study:

**At 4 weeks:** All samples were studied with a scanner; cement appeared homogeneous, with a high density but without any artifacts. The density size of the KING 0,01 cm<sup>2</sup> varied between 1680 and 1790 UH. The contours of the implants were clear and bone/cement contact appeared to be direct. (fig. 3).

The histology showed that the periphery of the cement was colonized by lately formed osseous spans, bordered by osteoplastic cells, which came to replace the cement gradually. Between the spans, the tissue was richly vascularized. There was no fibrous reaction between the newly formed bone and the mature bone located at the periphery of the cement. No inflammatory reaction was observed. (fig. 4).

At 11 weeks, no significant modification was noted with the scanner compared to the samples at 4 weeks. The cement always appeared dense (the measurements of density carried out in the center of the cement found values close to those of the first time of retreat). The contact with the bone was direct without edging of repair. First osteocondensation areas are present in the seat juxtamedullary cement.

#### Conclusion

The experimental plans were used to optimize the formulation of two types of injectable cement. In one case, it is necessary to introduce 1 % of silicone in weight into the solid phase (0.7 % in total weight, liquid plus solid phases) and in the other case 2 % of polyglycol in weight in the liquid phase (0.8 %) in total weight). The main physical chemical properties (evolution, mechanical properties, setting time) of the cements are preserved. The role of each additive is however slightly different. Silicone improves the rheology of cement, due to its lubricating character, but also slows down the rate of the evolution of the cement. Polyglycol has an effect on the rheology of cement without modifying its rate of the evolution. The first histological observations do not show fibro-inflammatory reactions. A large quantity of new formed bone appeared from the first month and osteointegration continued till the third month.

#### Acknowledgements

This work was completed in collaboration with the company Teknimed S.A established to 65500 Vic in Bigorre and supported by a CEE grant n° BMH4-97-9119.

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### Experimental increase of osteoporotic distal radius strength using a calcium phosphate cement

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**Abstract**

Diagnosis and treatment of osteoporosis are focused on demineralisation but bone mineral density is not directly correlated with bone strength. As with every material, the mechanical strength of bone depends upon its Young's modulus and its cross-sectional moment of inertia. In the clinical situation, bone strength can be quantified using peripheral quantitative computed tomography imaging (pQCT), a non-invasive imaging method, which allows calculation of a strength index. In this study, we tried to increase the fracture threshold of the distal radius by directly increasing bone strength rather than density. Twenty wrists in 10 cadavers were filled percutaneously with a calcium phosphate cement. Fluoroscopy and pQCT were performed twice, once before cementing and again 24 h after cement crystallisation to hydroxyapatite. We obtained measurements of trabecular and total bone density, and also stress strain index (SSI). Our results showed that trabecular bone density increased by a factor of 2.85, whereas total bone density increased by 1.61 and SSI by 1.99. Fluoroscopy showed two small leaks of cement at the point of injection. This study demonstrated that percutaneous injection of calcium phosphate cement increased distal radius strength, and consequently its fracture threshold. This technique could be employed in the future to prevent the occurrence of fractures in osteoporotic patients.

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Société General Electrics, Buc, France.

## Références

- [1] Riggs BL, Melton LJ. The prevention and treatment of osteoporosis. *New England Journal of Medicine* 1992;327:620-7.
- [2] Consensus development conference: diagnosis, prophylaxis and treatment of osteoporosis. *American Journal of Medicine*. 94. 1993. p. 646-50.
- [3] Nijs J, Westhovens R, Joly J, Cheng XG, Borghs H, Dequeker J. Diagnostic sensitivity of peripheral quantitative computed tomography measurements at ultra distal and proximal radius in postmenopausal women. *Bone* 1998;22(6):659-64.
- [4] Leichter I, Margulies JY, Weinreb A, et al. The relationship between bone density mineral content, and mechanical strength in the femoral neck. *Clinical Orthopaedics* 1982;163:272-81.
- [5] Liberman UA, Weiss SR, Bröll J, et al., The alendronate phase III osteoporosis treatment study group. Effect of oral alendronate on bone mineral density and the incidence of fractures in postmenopausal osteoporosis. *New England Journal of Medicine* 1995;333(22):1437-43.
- [6] Burger H, Van Daele PLA, Agra D, et al. Vertebral deformities as predictors of non-vertebral fractures. *British Medical Journal* 1994; 309:991-2.
- [7] Mc Clung MR. Clinical risk factors and evaluation of the risk of osteoporosis in clinical practice. *Ann Med Interne* 2000;151(5):392-8.
- [8] Schiessl H, Ferretti JL, Tysarczik-Niemeyer G, Willnecker J. Non-invasive bone strength index as analysed by peripheral quantitative computed tomography (pQCT). In: Schönau E, editor. *Paediatric osteology: new developments in diagnostics and therapy*. Elsevier Science; 1996. p. 141-6.
- [9] Ferretti JL, Capozza RF, Tysarczik-Niemeyer G, Schiessl H, Steffens M. Tomographic determination of stability parameters allow non-invasive estimation of bending or torsion strength. *Osteoporos Int* 1995;5:298-304.
- [10] Perren SM, Cordey J, Moncade N. Éléments de biomécanique de l'os et de l'ostéosynthèse. In: Kenesi C, editor. *Biomécanique générale*. Expansion Scientifique Française; 1985. p. 93-104.
- [11] Schoenau E, Neu C, Rauch F, Manz F. The development of bone strength at the proximal radius during childhood and adolescence. *The Journal of Clinical Endocrinology and Metabolism* 2001;86(2):613-8.
- [12] Ferretti JL, Capozza RF, Zanchetta JR. Mechanical validation of a tomographic (pQCT) index for non-invasive estimation of rat femur bending strength. *Bone* 1996;18(2):97-102.
- [13] Augat P, Reeb H, Claes LE. Prediction of fracture load at different skeletal sites by geometric properties of the cortical shell. *Journal of Bone and Mineral Research* 1996;11:1356-63.
- [14] Schneider PF, Fischer M, Allolio B, Felsenberg D, Schröder V, Semler J, Iitner J. Alendronate increases bone density and bone strength at the distal radius in postmenopausal women. *Journal of Bone and Mineral Research* 1999;14(8):1387-93.
- [15] Eastell R. Treatment of postmenopausal osteoporosis. *New England Journal of Medicine* 1998;338:736-46.
- [16] Black DM, Cummings SR, Karpf DB, et al., Fracture Intervention Trial Research Group. Randomised trial of effect of alendronate on risk of fracture in women with existing vertebral fractures. *The Lancet* 1996;348(9041):1535-41.

**Special 510(k)**

**Declaration of Conformity with Design Controls**

To the best of my knowledge, the verification activities, as required by the risk analysis, for the modification were performed by the designated individual(s) and the results demonstrated that the predetermined acceptance criteria were met.

_____	Date	_____ Name
<Title> P.R.D		GONCALVES Stephane
<Company> TEKNIUM	04.10.04	

The manufacturing facility, [Company Name] is in conformance with the design control requirements as specified in 21 CFR 820. 30 and the records are available for review.

_____	Date	_____ Name
<Title> P.R.D		GONCALVES Stephane
<Company> TEKNIUM	04.10.04	

**Differences Between Cementek® and Cementek® LV**  
**Effect on Risk**

Characteristic	Cementek®	Cementek® LV	Risk	Is Risk Increased?
<p><i>Material(s), including all additives</i></p>	<p>(b) (4)</p>			
<p><i>Physical form</i></p>				

Characteristic	Cementek®	Cementek® LV	Risk	Is Risk Increased?
<i>Mass, volume and density</i>	(b) (4)			
<i>Liquid: solid ratio</i>				

<b>Is Risk Increased?</b>	(b) (4)		
<b>Risk</b>	(b) (4)		
<b>Cementek® LV</b>	(b) (4)		
<b>Cementek®</b>	(b) (4)		
<b>Characteristic</b>	introduction to surgical site		

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
Dear Sir/Madam,

This supplement is in response to the deficiency letter received November 15, 2004.

1. A revised Design Control Activities Summary table is attached. In it all possible risks associated with each modification have been identified along with the verifications activities, the acceptance criteria and results.
2. A revised Information for Use leaflet is attached. It addresses the concerns regarding adding additional substances, excess material, contraindication for vertebral compression fractures and over-pressurizing the material.

I trust that these deficiencies have been adequately addressed.

Respectfully,

J.D. Webb  
Authorized Contact Person

65

**Design Control Activities Summary**  
**Cementek® LV**

<b>Modification</b>	<b>Risk</b>	<b>Verification Activity</b>	<b>Acceptance Criteria</b>	<b>Results of Verification</b>
(b) (4)				

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Cementek LV® is contraindicated for treating vertebral compression fractures.

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Do not add any additional substances other than those provided in this package.

Avoid over-pressurizing Cementek LV® as this may lead to extrusion of the device beyond the site of its intended applications and damage the surrounding tissues.

Over pressurizing may also lead to fat embolization or embolization of the material into the bloodstream.

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**POSSIBLE ADVERSE EFFECTS**

A successful result is not achieved in every surgical case. Re-operation to remove or replace the implant may be required

67

at any time due to medical reasons or device failure. One or more of the following complications may occur and corrective action should be taken:

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**STERILIZATION**

Cementek ® LV is sterilized with 25 to 40 kGy (2,5 to 4,0 Mrad) gamma radiation. All sterile implants are to be kept unopened in the original packaging until the time of implantation. Before using the implant, the protective packaging should be checked for damage as this could be detrimental to the sterility. Aseptic procedures should be observed when removing the implant from the protective package.

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**Remarks:**

Verify the integrity of the packaging before use, the guarantee of sterilization is 5 years from the date of sterilization. Use after the peremptory date is not allowed. All re-sterilization of the product is forbidden use only once.

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Public Health Service  
Food and Drug Administration  
Memorandum

From: Reviewer(s) - Name(s) *Valerie J. Green*

Subject: 510(k) Number *K042911/51*

To: The Record - It is my recommendation that the subject 510(k) Notification:

- Refused to accept.
- Requires additional information (other than refuse to accept).
- Is substantially equivalent to marketed devices.
- NOT substantially equivalent to marketed devices.
- Other (e.g., exempt by regulation, not a device, duplicate, etc.).

- Is this device subject to Section 522 Postmarket Surveillance?  YES  NO
- Is this device subject to the Tracking Regulation?  YES  NO
- Was clinical data necessary to support the review of this 510(k)?  YES  NO
- Is this a prescription device?  YES  NO
- Was this 510(k) reviewed by a Third Party?  YES  NO
- Special 510(k)?  YES  NO
- Abbreviated 510(k)? Please fill out form on H Drive 510k/boilers  YES  NO

- Truthful and Accurate Statement  Requested  Enclosed
- A 510(k) summary OR  A 510(k) statement
- The required certification and summary for class III devices *NA*
- The indication for use form

Combination Product Category (Please see algorithm on H drive 510k/Boilers) N

Animal Tissue Source  YES  NO Material of Biological Origin  YES  NO

The submitter requests under 21 CFR 807.95 (doesn't apply for SEs):  
 No Confidentiality  Confidentiality for 90 days  Continued Confidentiality exceeding 90 days

Predicate Product Code with class: \_\_\_\_\_ Additional Product Code(s) with panel (optional): \_\_\_\_\_

Review: *[Signature]*  
(Branch Chief)

*REDB*  
(Branch Code)

*12/23/04*  
(Date)

Final Review: *[Signature]*  
(Division Director)

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Public Health Service  
Food and Drug Administration  
Memorandum

From: Reviewer(s) - Name(s) *Nadine Sloan*  
Subject: 510(k) Number *K042911*

To: The Record - It is my recommendation that the subject 510(k) Notification:

- Refused to accept.
- Requires additional information (other than refuse to accept). *11/15/04*
- Is substantially equivalent to marketed devices.
- NOT substantially equivalent to marketed devices.
- Other (e.g., exempt by regulation, not a device, duplicate, etc.)

- |   |                              |                             |
|---|------------------------------|-----------------------------|
| Is this device subject to Section 522 Postmarket Surveillance?    | <input type="checkbox"/> YES | <input type="checkbox"/> NO |
| Is this device subject to the Tracking Regulation?                | <input type="checkbox"/> YES | <input type="checkbox"/> NO |
| Was clinical data necessary to support the review of this 510(k)? | <input type="checkbox"/> YES | <input type="checkbox"/> NO |
| Is this a prescription device?                                    | <input type="checkbox"/> YES | <input type="checkbox"/> NO |
| Was this 510(k) reviewed by a Third Party?                        | <input type="checkbox"/> YES | <input type="checkbox"/> NO |
| Special 510(k)?   | <input type="checkbox"/> YES | <input type="checkbox"/> NO |
| Abbreviated 510(k)? Please fill out form on H Drive 510k/boilers  | <input type="checkbox"/> YES | <input type="checkbox"/> NO |

- Truthful and Accurate Statement  Requested  Enclosed
- A 510(k) summary OR  A 510(k) statement
- The required certification and summary for class III devices
- The indication for use form

Combination Product Category (Please see algorithm on H drive 510k/Boilers) *N*

Animal Tissue Source  YES  NO Material of Biological Origin  YES  NO

The submitter requests under 21 CFR 807.95 (doesn't apply for SEs):  
 No Confidentiality  Confidentiality for 90 days  Continued Confidentiality exceeding 90 days

Predicate Product Code with class: \_\_\_\_\_ Additional Product Code(s) with panel (optional): \_\_\_\_\_

Review: *[Signature]* *R&OB* *11/15/04*  
(Branch Chief) (Branch Code) (Date)

Final Review: \_\_\_\_\_  
(Division Director) (Date)

Revised:4/2/03

**REVIEW MEMORANDUM**

**DATE:** 11/15/04

**FROM:** Nadine Y. Sloan  
DGRND/REDB

**SUBJECT:** K042911; Teknimed; Cementek LV Bone Substitute – Special 510(k)

**TO:** File

**Contact:** J.D. Webb 512-338-0199

---

This 510(k) is submitted as a special 510(k) requesting changes to the company's previously cleared bone void filler (Cementek Bone Substitute; K041493).

(b) (4)

Recommended Deficiencies:

(b) (4)

(b) (4)


  
Madeline Y. Sloan  
cementek.Special.AI.rev

## Internal Administrative Form

	YES	NO
1. Did the firm request expedited review?		
2. Did we grant expedited review?		
3. Have you verified that the Document is labeled Class III for GMP purposes?	NA	
4. If, not, has POS been notified?		
5. Is the product a device?	✓	
6. Is the device exempt from 510(k) by regulation or policy?	✓	✓
7. Is the device subject to review by CDRH?		✓
8. Are you aware that this device has been the subject of a previous NSE decision?		✓
9. If yes, does this new 510(k) address the NSE issue(s), (e.g., performance data)?		
10. Are you aware of the submitter being the subject of an integrity investigation?		✓
11. If, yes, consult the ODE Integrity Officer.		
12. Has the ODE Integrity Officer given permission to proceed with the review? (Blue Book Memo #191-2 and Federal Register 90N0332, September 10, 1991.		

The deficiencies identified above represent the issues that we believe need to be resolved before our review of your 510(k) submission can be successfully completed. In developing the deficiencies, we carefully considered the statutory criteria as defined in Section 513(i) of the Federal Food, Drug, and Cosmetic Act for determining substantial equivalence of your device. We also considered the burden that may be incurred in your attempt to respond to the deficiencies. We believe that we have considered the least burdensome approach to resolving these issues. If, however, you believe that information is being requested that is not relevant to the regulatory decision or that there is a less burdensome way to resolve the issues, you should follow the procedures outlined in the "A Suggested Approach to Resolving Least Burdensome Issues" document. It is available on our Center web page at: <http://www.fda.gov/cdrh/modact/leastburdensome.html>

### 510(k) "SUBSTANTIAL EQUIVALENCE" DECISION-MAKING PROCESS

\* 510(k) Submissions compare new devices to marketed devices. FDA requests additional information if the relationship between marketed and "predicate" (pre-Amendments or reclassified post-Amendments) devices is unclear.

\*\* This decision is normally based on descriptive information alone, but limited testing information is sometimes required.

\*\*\* Data maybe in the 510(k), other 510(k)s, the Center's classification files, or the literature.

## Internal Administrative Form

	YES	NO
1. Did the firm request expedited review?	<i>Special</i> ✓	
2. Did we grant expedited review?	NA	
3. Have you verified that the Document is labeled Class III for GMP purposes?	✓	
4. If, not, has POS been notified?	✓	✓
5. Is the product a device?	✓	
6. Is the device exempt from 510(k) by regulation or policy?	✓	
7. Is the device subject to review by CDRH?		✓
8. Are you aware that this device has been the subject of a previous NSE decision?		✓
9. If yes, does this new 510(k) address the NSE issue(s), (e.g., performance data)?		✓
10. Are you aware of the submitter being the subject of an integrity investigation?		✓
11. If, yes, consult the ODE Integrity Officer.		
12. Has the ODE Integrity Officer given permission to proceed with the review? (Blue Book Memo #191-2 and Federal Register 90N0332, September 10, 1991.		

**SPECIAL 510(k): Device Modification**  
**ODE Review Memorandum (Decision Making Document is Attached)**

**To:** THE FILE

**RE:** DOCUMENT NUMBER K042911  
Cementek LV (Teknimed, S.A.)

---

This 510(k) submission contains information/data on modifications made to the SUBMITTER'S own Class II, Class III or Class I devices requiring 510(k). The following items are present and acceptable (delete/add items as necessary):

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510(k) Number: K042911  
Product: CEMENTEK LV

The additional information you have submitted has been received.

We will notify you when the processing of this submission has been completed or if any additional information is required. Please remember that all correspondence concerning your submission MUST be sent to the Document Mail Center (HFZ-401) at the above letterhead address. Correspondence sent to any address other than the one above will not be considered as part of your official premarket notification submission. Also, please note the new Blue Book Memorandum regarding Fax and E-mail Policy entitled, "Fax and E-Mail Communication with Industry about Premarket Files Under Review. Please refer to this guidance for information on current fax and e-mail practices at [www.fda.gov/cdrh/ode/a02-01.html](http://www.fda.gov/cdrh/ode/a02-01.html).

The Safe Medical Devices Act of 1990, signed on November 28, states that you may not place this device into commercial distribution until you receive a letter from FDA allowing you to do so. As in the past, we intend to complete our review as quickly as possible. Generally we do so 90 days. However, the complexity of a submission or a requirement for additional information may occasionally cause the review to extend beyond 90 days. Thus, if you have not received a written decision or been contacted within 90 days of our receipt date you may want to check with FDA to determine the status of your submission.

If you have procedural or policy questions, please contact the Division of Small Manufacturers International and Consumer Assistance (DSMICA) at (301) 443-6597 or at their toll-free number (800) 638-2041, or contact me at (301) 594-1190.

Sincerely yours,

Marjorie Shulman  
Supervisory Consumer Safety Officer  
Premarket Notification Section  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

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
Dear Sir/Madam,

This supplement is in response to the deficiency letter received November 15, 2004.

1. A revised Design Control Activities Summary table is attached. In it all possible risks associated with each modification have been identified along with the verifications activities, the acceptance criteria and results.
2. A revised Information for Use leaflet is attached. It addresses the concerns regarding adding additional substances, excess material, contraindication for vertebral compression fractures and over-pressurizing the material.

I trust that these deficiencies have been adequately addressed.

Respectfully,

J.D. Webb  
Authorized Contact Person

*5132*

*12*

**Design Control Activities Summary**  
**Cementek® LV**

<b>Modification</b>	<b>Risk</b>	<b>Verification Activity</b>	<b>Acceptance Criteria</b>	<b>Results of Verification</b>
(b) (4)	(b) (4)	(b) (4)	(b) (4)	(b) (4)

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Cementek LV® is contraindicated for treating vertebral compression fractures.

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Do not add any additional substances other than those provided in this package.

Avoid over-pressurizing Cementek LV® as this may lead to extrusion of the device beyond the site of its intended applications and damage the surrounding tissues.

Over pressurizing may also lead to fat embolization or embolization of the material into the bloodstream.

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**POSSIBLE ADVERSE EFFECTS**

A successful result is not achieved in every surgical case. Re-operation to remove or replace the implant may be required

at any time due to medical reasons or device failure. One or more of the following complications may occur and corrective action should be taken:

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#### **STERILIZATION**

Cementek ® LV is sterilized with 25 to 40 kCiy (2,5 to 4,0 Mrad) gamma radiation. All sterile implants are to be kept unopened in the original packaging until the time of implantation. Before using the implant, the protective packaging should be checked for damage as this could be detrimental to the sterility. Aseptic procedures should be observed when removing the implant from the protective package.

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Store at room temperature

**Generic "Design Control Activities Summary"**

Device Modification	Risk	Verification Activity	Acceptance Criteria	Results of Verification
---------------------	------	-----------------------	---------------------	-------------------------

(b) (4)

**Sample "Design Control Activities Summary"**

Modification	Risk	Verification Activity	Acceptance Criteria	Results of Verification
--------------	------	-----------------------	---------------------	-------------------------

(b) (4)

16

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© 2003 Elsevier SAS. Tous droits réservés.