



U.S. Department of Health & Human Services

Food and Drug Administration

SAVE REQUEST

USER: (kml)
FOLDER: K063166 - 217 pages
COMPANY: I.T.S. IMPLANTAT-TECHNOLOGIE-SYSTEME GMBH (ITSIMPL)
PRODUCT: PLATE, FIXATION, BONE (HRS)
SUMMARY: Product: I.T.S. PELVIC RECONSTRUCTION SYSTEM

DATE REQUESTED: Apr 27, 2016

DATE PRINTED: Apr 27, 2016

Note: Printed



K063166

Section XII: 510(k) Summary of Safety and Effectiveness

SAFE MEDICAL DEVICES ACT OF 1990
510(k) Summary

DEC 22 2006

NAME OF FIRM: I.T.S. Implantat-Technologie-Systeme GmbH.
Autal 28.
Lassnitzhoehe A – 8301
AUSTRIA

510(k) FIRM CONTACT: Al Lippincott
Engineering Consulting Services, Inc.
3150 E. 200th St.
Prior Lake, MN 55372

TRADE NAME: Pelvic Reconstruction System (PRS)

COMMON NAME: Pelvic Bone Fixation Set

CLASSIFICATION: Plate, Fixation, Bone (see 21 CFR, Sec. 888.3030),
Pin, Fixation, Threaded (see 21 CFR, Sec. 888.3040)
Screw, Fixation, Bone (see 21 CFR, Sec. 888.3040),
Washer, Bolt, Nut (see 21 CFR, Sec. 888.3030)

DEVICE PRODUCT CODE: HRS

SUBSEQUENT PRODUCT CODE: JDW, HWC, HTN

SUBSTANTIALLY EQUIVALENT DEVICES Synthes 3.5mm Low Profile Pelvic Reconstruction (**K031573**)
Stryker Trauma Pelvic Set (**K001614**)
Zimmer Pelvic Reconstruction Set
Synthes Sacral Bar System (**K001720**)
I.T.S. GmbH FR.O.H. Calcaneus Repair System (**K051642**)
Synthes 3.9mm Pelvic Screw (**K013044**)
I.T.S. Straight Plate with Angular Stability & Screw System (**K060156**)

DEVICE DESCRIPTION: The I.T.S. Pelvic Reconstruction System (PRS) encompasses a number of fracture fixation subsystems (multiple pelvic plate designs, sacral threaded rod, and cannulated screw & washer) for fracture fixation and reconstruction of pelvic ring fractures in the pelvis.

The I.T.S. PRS Low Profile – Multiple Type – Pelvic Plating System consists of the following plate types: 1) A Straight Plate at a 2.5mm thickness with 10 to 14 hole length sizes, 2) A Curved Plate with a 108mm radius at a 2.5mm thickness with 4 to 16 hole length sizes, 3) A Curved Plate with a 88mm radius at a 2.5mm thickness also with 4 to 16 hole length sizes, 4) A Symphysis Plate at a 4.0mm thickness in both a 4 & 6 hole size, 5) A Sacroiliac-Joint (SIJ) L-Shaped Plate at a 2.5mm thickness in a left and right 5 hole size, and 6) A SIJ Closed Plate at a 2.5mm thickness in a 4 hole size.

All plate designs are low profile in thickness and made from CP Titanium. The PRS Pelvic Plating System also encompasses a number

of Cancellous (5.9mm std.compression & 5.9mm Locking screws) and Cortical (4.5mm std.compression screw) screw types and length sizes. All bone screws are pre-drilling and self-tapping in design and manufactured from high strength 6-4 Alloyed Titanium. All components (plates & screws) have a TIODIZE II surface treatment preparation.

The I.T.S. PRS Sacral Rod System consists of a threaded pin, wedge-shaped washers, and nut/locknut design with one end of the threaded pin having a trocar point. The assembled unit is used to stabilize ilio-iliac posterior pelvic ring disruption injuries. All components are manufactured from high strength 6-4 Alloyed Titanium material and have the TIODIZE II surface preparation.

The I.T.S. PRS 7.3mm Cannulated SIJ Traction Screw & Washer is a pre-drilling, self-tapping, and back-tapping screw design for guided reduction of pelvic bone fractures. Washers (flat & curved) are available for use with the screw. All 7.3mm Cannulated Screws & washers are manufactured from high strength 6-4 Alloyed Titanium and have the TIODIZE II surface preparation.

INTENDED USE:

The intended use of the I.T.S. Pelvic Reconstruction System (PRS) is to stabilize one or more pelvic bone fractures in the pelvic ring area of the pelvis of an adult or pediatric patient which include the use of the following pelvic fracture fixation systems.

Indications for use of the I.T.S. PRS Low Profile - Multiple Type – Pelvic Plating System include:

- 1) Fracture reconstruction of the acetabulum, sacrum, ilium, and entire pelvic ring,
- 2) Revision surgery of pseudarthroses, non-unions and mal-unions,
- 3) Ilio-Iliac distance osteosynthesis,
- 4) Osteotomies,
- 5) Arthrodesis,
- 6) Sacroiliac joint dislocations, and
- 7) Symphysis pubis disruptions

Indications for use of the I.T.S. PRS Sacral Rod System include:

- 1) Fixation of fractures of the posterior pelvis,
- 2) Fixation of fractures of the posterior iliac spine,
- 3) Fixation of fractures of the posterior inferior iliac spine,
- 4) Dorsal stabilization of the posterior pelvic ring for unstable pelvic ring injuries,
- 5) Fixation of sacral fractures, and
- 6) Fracture dislocations of the sacro-iliac joint

Indications for use of the I.T.S. PRS 7.3mm Cannulated SIJ Traction Screw & Washer include fracture fixation of pelvic bone fractures where indicated and for periacetabular osteotomies.

The system(s) has not been studied in spinal use, and is not intended for use in vertebral column fracture or fusion procedures.

BASIS OF SUBSTANTIAL EQUIVALENCE:

The I.T.S. Pelvic Reconstruction System (PRS) is substantially equivalent to the Synthes, Zimmer, and Stryker pelvic reconstruction systems and the I.T.S. GmbH stabilizing bone plating systems.

SUMMARY OF SAFETY AND EFFECTIVENESS:

The I.T.S. Pelvic Reconstruction System (PRS) is shown to be safe and effective for use in 'pelvic ring' bone fracture fixation of the pelvis.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

I.T.S Implantat-Technologie-Systeme GmbH
% Mr. Al Lippincott
Engineering Consulting Services, Inc.
3150 E. 200th Street
Prior Lake, Minnesota 55372

DEC 22 2006

Re: K063166

Trade/Device Name: Pelvic Reconstruction System (PRS)

Regulation Number: 21 CFR 888.3030

Regulation Name: Single/multiple component metallic bone fixation appliances and accessories

Regulatory Class: Class II

Product Code: HRS, JDW, HWC

Dated: October 13, 2006

Received: October 18, 2006

Dear Mr. Lippincott:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Mr. Al Lippincott

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0210. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson". The signature is written in a cursive style with a small "to" written below the main signature.

Mark N. Melkerson
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure



Indications for Use

510(k) NUMBER: K063166

DEVICE NAME: PELVIC RECONSTRUCTION SYSTEM (PRS)

The ***intended use*** of the I.T.S. Pelvic Reconstruction System (PRS) is to stabilize one or more pelvic bone fractures in the 'pelvic ring' area of the pelvis of an adult or pediatric patient which include the use of the following pelvic fracture fixation systems.

Indications for use of the I.T.S. PRS Low Profile - Multiple Type - Pelvic Plating System include:

- 1) Fracture reconstruction of the acetabulum, sacrum, ilium, and entire pelvic ring,
- 2) Revision surgery of pseudarthroses, non-unions and mal-unions,
- 3) Ilio-Iliac distance osteosynthesis,
- 4) Osteotomies,
- 5) Arthrodesis,
- 6) Sacroiliac joint dislocations, and
- 7) Symphysis pubis disruptions

Indications for use of the I.T.S. PRS Sacral Rod System include:

- 1) Fixation of fractures of the posterior pelvis,
- 2) Fixation of fractures of the posterior iliac spine,
- 3) Fixation of fractures of the posterior inferior iliac spine,
- 4) Dorsal stabilization of the posterior pelvic ring for unstable pelvic ring injuries,
- 5) Fixation of sacral fractures, and
- 6) Fracture dislocations of the sacro-iliac joint

Indications for use of the I.T.S. PRS 7.3mm Cannulated SIJ Traction Screw & Washer include fracture fixation of pelvic bone fractures where indicated and for periacetabular osteotomies.

The system(s) has not been studied in spinal use, and is not intended for use in vertebral column fracture or fusion procedures.

Prescription Use X AND/OR Over-The-Counter-Use _____
 (Per 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Barbara Buchard
 (Division Sign-Off)

Division of General, Restorative,
 and Neurological Devices

Concurrence of CDRH, Office of Device Evaluation (ODE)

510(k) Number K063166



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

I.T.S Implantat-Technologie-Systeme GmbH
% Mr. Al Lippincott
Engineering Consulting Services, Inc.
3150 E. 200th Street
Prior Lake, Minnesota 55372

DEC 22 2006

Re: K063166

Trade/Device Name: Pelvic Reconstruction System (PRS)
Regulation Number: 21 CFR 888.3030
Regulation Name: Single/multiple component metallic bone fixation appliances and accessories
Regulatory Class: Class II
Product Code: HRS, JDW, HWC
Dated: October 13, 2006
Received: October 18, 2006

Dear Mr. Lippincott:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Mr. Al Lippincott

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0210 Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson". The signature is written in a cursive style with a large initial "M".

Mark N. Melkerson
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) NUMBER: K063166

DEVICE NAME: PELVIC RECONSTRUCTION SYSTEM (PRS)

The **intended use** of the I.T.S. Pelvic Reconstruction System (PRS) is to stabilize one or more pelvic bone fractures in the 'pelvic ring' area of the pelvis of an adult or pediatric patient which include the use of the following pelvic fracture fixation systems.

Indications for use of the I.T.S. PRS Low Profile - Multiple Type - Pelvic Plating System include:

- 1) Fracture reconstruction of the acetabulum, sacrum, ilium, and entire pelvic ring,
- 2) Revision surgery of pseudarthroses, non-unions and mal-unions,
- 3) Ilio-Iliac distance osteosynthesis,
- 4) Osteotomies,
- 5) Arthrodesis,
- 6) Sacroiliac joint dislocations, and
- 7) Symphysis pubis disruptions

Indications for use of the I.T.S. PRS Sacral Rod System include:

- 1) Fixation of fractures of the posterior pelvis,
- 2) Fixation of fractures of the posterior iliac spine,
- 3) Fixation of fractures of the posterior inferior iliac spine,
- 4) Dorsal stabilization of the posterior pelvic ring for unstable pelvic ring injuries,
- 5) Fixation of sacral fractures, and
- 6) Fracture dislocations of the sacro-iliac joint

Indications for use of the I.T.S. PRS 7.3mm Cannulated SIJ Traction Screw & Washer include fracture fixation of pelvic bone fractures where indicated and for periacetabular osteotomies.

The system(s) has not been studied in spinal use, and is not intended for use in vertebral column fracture or fusion procedures.

Prescription Use X AND/OR Over-The-Counter-Use _____

(Per 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Barbara Pruchno
(Division Sign-Off)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Division of General, Restorative,
and Neurological Devices

510(k) Number K063166

Food and Drug Administration
Center for Devices and
Radiological Health
Office of Device Evaluation
Document Mail Center (HFZ-401)
9200 Corporate Blvd.
Rockville, Maryland 20850

October 18, 2006

I.T.S. IMPLANTAT-TECHNOLOGIE-SYSTEM 510(k) Number: K063166
C/O ENGINEERING CONSULTING SERVICES Received: 18-OCT-2006
3150 E. 200TH ST. Product: I.T.S. PELVIC
PRIOR LAKE, MN 55372 RECONSTRUCTION
ATTN: AL LIPPINCOTT SYSTEM

The Food and Drug Administration (FDA), Center for Devices and Radiological Health (CDRH), has received the Premarket Notification, (510(k)), you submitted in accordance with Section 510(k) of the Federal Food, Drug, and Cosmetic Act (Act) for the above referenced product and for the above referenced 510(k) submitter. Please note, if the 510(k) submitter is incorrect, please notify the 510(k) Staff immediately. We have assigned your submission a unique 510(k) number that is cited above. Please refer prominently to this 510(k) number in all future correspondence that relates to this submission. We will notify you when the processing of your 510(k) has been completed or if any additional information is required. YOU MAY NOT PLACE THIS DEVICE INTO COMMERCIAL DISTRIBUTION UNTIL YOU RECEIVE A LETTER FROM FDA ALLOWING YOU TO DO SO.

Please remember that all correspondence concerning your submission MUST be sent to the Document Mail Center (DMC) (HFZ-401) at the above letterhead address. Correspondence sent to any address other than the one above will not be considered as part of your official 510(k) submission.

Please note the following documents as they relate to 510(k) review:
1) Guidance for Industry and FDA Staff entitled, "FDA and Industry Actions on Premarket Notification (510(k)) Submissions: Effect on FDA Review Clock and Performance Assessment". The purpose of this document is to assist agency staff and the device industry in understanding how various FDA and industry actions that may be taken on 510(k)s should affect the review clock for purposes of meeting the Medical Device User Fee and Modernization Act (MDUFMA). Please review this document at www.fda.gov/cdrh/mdufma/guidance/1219.html.
2) Guidance for Industry and FDA Staff entitled, "Format for Traditional and Abbreviated 510(k)s". This guidance can be found at www.fda.gov/cdrh/ode/guidance/1567.html. Please refer to this guidance for assistance on how to format an original submission for a Traditional or Abbreviated 510(k).
3) Blue Book Memorandum regarding Fax and E-mail Policy entitled, "Fax and E-Mail Communication with Industry about Premarket Files Under Review". Please refer to this guidance for information on current fax and e-mail practices at www.fda.gov/cdrh/ode/a02-01.html.

In all future premarket submissions, we encourage you to provide an electronic copy of your submission. By doing so, you will save FDA resources and may help reviewers navigate through longer documents more easily. Under CDRHs e-Copy Program, you may replace one paper copy of an premarket submission (e.g., 510(k), IDE, PMA, HDE) with an electronic copy. For more information about the program, including the formatting requirements, please visit our web site at www.fda.gov/cdrh/elecsb.html.

Lastly, you should be familiar with the regulatory requirements for medical devices available at Device Advice www.fda.gov/cdrh/devadvice/. If you have questions on the status of your submission, please contact DSMICA at (240) 276-3150 or the toll-free number (800) 638-2041, or at their Internet address <http://www.fda.gov/cdrh/dsma/dsmastaf.html>. If you have policy or procedural questions, please contact anyone on the 510(k) Staff at (240)276-4040.

Sincerely yours,

Marjorie Shulman
Supervisory Consumer Safety Officer
Office of Device Evaluation
Center for Devices and Radiological Health

K063 166

Engineering Consulting Services, Inc.

3150 E. 200th St.
Prior Lake, MN. 55372

(952) 492-5858 TEL
(952) 492-5859 FAX

[Handwritten signature]

October 13, 2006

Food and Drug Administration
Center for Devices and Radiological Health
Document Mail Center, HFZ-401
1390 Piccard Drive
Rockville, MD 20850

K-26

Re: 510(k) Notification – I.T.S. GmbH Pelvic Reconstruction System (PRS)

Dear Sir/Madam,

Enclosed is a submission pursuant to Section 510(k) of the Federal Food, Drug, and Cosmetic Act and the regulations contained in Title 21 CFR 807. This submission is to notify the FDA of our intent to market the I.T.S. Pelvic Reconstruction System (PRS). We believe this fracture fixation system to be substantially equivalent to other pelvic bone fixation devices that are commercially available in the U.S..

Attached is the FDA Form # 3601 and copy of Check # [redacted] showing the 'User Fee' payment with PIN [redacted] via overnight Registered Priority Mail on 10-10-2006.

Some of the information contained herein is considered **PROPRIETARY** and is stamped accordingly.

Included in this submission in Section I is a Truthful and Accurate Statement, in Section XI is an Indications Enclosure, and in Section XII is a 510(K) Summary.

Please direct any inquiries concerning this submission to **Al Lippincott**, Engineering Consulting Services, Inc., U.S. Agent and Official Correspondent to I.T.S. Implantat-Technologie-Systeme GmbH, at **(952) 492-5858** or FAX **(952) 492-5859**.

Thank you for your prompt attention to this submission.

Sincerely,

Al Lippincott

Al Lippincott
U.S. Agent and Official Correspondent for I.T.S. Implantat-Technologie-Systeme GmbH

Enclosure: 510(k) Submission - I.T.S. GmbH - Pelvic Reconstruction System (PRS)

K26
OR
II

16

DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION MEDICAL DEVICE USER FEE COVER SHEET	PAYMENT IDENTIFICATION NUMBER: b(4) Write the Payment Identification number on your check.
---	---

A completed Cover Sheet must accompany each original application or supplement subject to fees. The following actions must be taken to properly submit your application and fee payment:

1. Electronically submits the completed Cover Sheet to the Food and Drug Administration (FDA) before payment is sent.
2. Include printed copy of this completed Cover Sheet with a check made payable to the Food and Drug Administration. Remember that the Payment Identification Number must be written on the check.
3. Mail Check and Cover Sheet to the US Bank Lock Box, FDA Account, P.O. Box 956733, St. Louis, MO 63195-6733. (Note: In no case should payment be submitted with the application.)
4. If you prefer to send a check by a courier, the courier may deliver the check and Cover Sheet to: US Bank, Attn: Government Lockbox 956733, 1005 Convention Plaza, St. Louis, MO 63101. (Note: This address is for courier delivery only. Contact the US Bank at 314-418-4821 if you have any questions concerning courier delivery.)
5. For Wire Transfer Payment Procedures, please refer to the MDUFMA Fee Payment Instructions at the following URL: <http://www.fda.gov/cdrh/mdufma/faqs.html#3a>. You are responsible for paying all fees associated with wire transfer.
6. Include a copy of the complete Cover Sheet in volume one of the application when submitting to the FDA at either the CBER or CDRH Document Mail Center.

1. COMPANY NAME AND ADDRESS (include name, street address, city state, country, and post office code) ITS IMPLANT TECHNICAL SYSTEMS GMBH (ASUTRIA) 3150 e. 200th. St. Prior Lake MN 55372 US 1.1 EMPLOYER IDENTIFICATION NUMBER (EIN) 411825819	2. CONTACT NAME Albert Lippincott III 2.1 E-MAIL ADDRESS allippincott@msn.com 2.2 TELEPHONE NUMBER (include Area code) 952-492-5858 2.3 FACSIMILE (FAX) NUMBER (Include Area code) 952-492-5859
---	--

3. TYPE OF PREMARKET APPLICATION (Select one of the following in each column; if you are unsure, please refer to the application descriptions at the following web site: <http://www.fda.gov/dc/mdufma>)

Select an application type: <input checked="" type="checkbox"/> Premarket notification(510(k)); except for third party <input type="checkbox"/> Biologics License Application (BLA) <input type="checkbox"/> Premarket Approval Application (PMA) <input type="checkbox"/> Modular PMA <input type="checkbox"/> Product Development Protocol (PDP) <input type="checkbox"/> Premarket Report (PMR)	3.1 Select one of the types below <input checked="" type="checkbox"/> Original Application Supplement Types: <input type="checkbox"/> Efficacy (BLA) <input type="checkbox"/> Panel Track (PMA, PMR, PDP) <input type="checkbox"/> Real-Time (PMA, PMR, PDP) <input type="checkbox"/> 180-day (PMA, PMR, PDP)
--	---

4. ARE YOU A SMALL BUSINESS? (See the instructions for more information on determining this status)

YES, I meet the small business criteria and have submitted the required qualifying documents to FDA NO, I am not a small business

4.1 If Yes, please enter your Small Business Decision Number:

5. IS THIS PREMARKET APPLICATION COVERED BY ANY OF THE FOLLOWING USER FEE EXCEPTIONS? IF SO, CHECK THE APPLICABLE EXCEPTION.


<input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates, parents, and partner firms <input type="checkbox"/> This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only	<input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population <input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially
---	--

6. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION FOR WHICH FEES WERE WAIVED DUE TO SOLE USE IN A PEDIATRIC POPULATION THAT NOW PROPOSES CONDITION OF USE FOR ANY ADULT POPULATION? (If so, the application is subject to the fee that applies for an original premarket approval application (PMA).)

YES NO

7. USER FEE PAYMENT AMOUNT SUBMITTED FOR THIS PREMARKET APPLICATION (FOR FISCAL YEAR 2005)

b(4) 09-Oct-2006



17



ENGINEERING CONSULTING SERVICES, INC.

3150 E. 200TH STREET
PRIOR LAKE, MN 55372
(952) 492-5858

(b)(4)

(b)(4)

10/9/2006

PAY TO THE ORDER OF **Food and Drug Administration**

(b)(4)

\$

(b)(4)

ENGINEERING CONSULTING SERVICES, INC.

Food and Drug Administration
Consulting: FDA 'User Fee'

10/9/2006

510(k) Application Fee for PRS for I.T.S. Austria

(b)(4)

(b)(4)

Business Checking -

(b)(4)

(b)(4)

ENGINEERING CONSULTING SERVICES, INC.

Food and Drug Administration
Consulting: FDA 'User Fee'

10/9/2006

510(k) Application Fee for PRS for I.T.S. Austria

(b)(4)

(b)(4)

PAYMENT RECORD

Business Checking -

(b)(4)

(b)(4)

18

510(k) Submission
I.T.S. Implantat-Technologie-Systeme GmbH
Pelvic Reconstruction System (PRS)


Table of Contents:

ADMINISTRATIVE INFORMATION

- Section I Truthful and Accurate Statement
- Section II Manufacturer Identification
- Section III Device Identification
- Section IV Device Background

DEVICE INFORMATION

- Section V Intended Uses/Indications
- Section VI Device Description
 - Drawing(s) of the Device(s)
- Section VII Materials
 - ASTM Specifications F 67 & F 136
 - Sample Certifications of CP and 6-4 alloyed Titanium
- Section VIII Testing Information
 - Engineering Rationale for not testing
 - Synthes, Zimmer, and Stryker marketing and literature references
- Section IX Sterilization Information
 - AAMI Guidelines for Sterilization
- Section X Labeling
 - Pelvic Reconstuction System (PRS) sample labeling
- Section XI Indications Enclosure
- Section XII 510(K) Summary or Statement
- Section XIII Substantial Equivalence Information
 - Comparison Table
 - 510(k) Summary Statements of equivalent devices with literature
- Section XIV Establishment Registration and Device Listing




PREMARKET NOTIFICATION
TRUTHFUL AND ACCURATE STATEMENT

(As required By 21 CFR 807.87(j))

I certify, in my capacity as president of I.T.S. Implantat-Technologie-Systeme GmbH, I believe to the best of my knowledge, that all data and information submitted in the premarket notification are truthful and accurate and that no material fact has been omitted.

Signed


.....
Dr. Eva Ruprechter

06/10/2006 Oct. 6, 2006
.....
Dated

.....
Premarket Notification (510(k)) Number

510(k) Submission
I.T.S. Implantat-Technologie-Systeme GmbH
Pelvic Reconstruction System (PRS)

Section II: Manufacturer Identification

1. Name and Address of Manufacturer:

I.T.S. Implantat-Technologie-Systeme GmbH
Autal 28
Lassnitzhoehe A – 8301

AUSTRIA

2. Sponsor of the 510(k) submission:

I.T.S. Implantat-Technologie-Systeme GmbH
Autal 28
Lassnitzhoehe A – 8301

AUSTRIA

3. Manufacturer's Registration Number:

Registration No.: **3004369035**

Owner/Operator No.: **9059922**

4. Contact Person for All Correspondence:

Al Lippincott,

U.S. Agent and Official Correspondent for I.T.S. Implantat-Technologie-Systeme GmbH

Engineering Consulting Services, Inc.
3150 E. 200th. St.
Prior Lake, MN. 55372

Telephone Number: (952) 492-5858

FAX Number: (952) 492-5859

Cell Number: (b) (6)

E-mail: allippincott@msn.com

510(k) Submission
I.T.S. Implantat-Technologie-Systeme GmbH
Pelvic Reconstruction System (PRS)

Section III: Device Identification

1. **Proprietary Name:** I.T.S. Pelvic Reconstruction System (PRS)

2. **Common Name:** Pelvic Bone Fixation Set

3. **Classification Name and Reference:**

Plate, Fixation, Bone (see 21 CFR, Sec. 888.3030)

Pin, Fixation, Threaded (see 21 CFR, Sec. 888.3040)

Screw, Fixation, Bone (see 21CFR, Sec. 888.3040)

Washer, Bolt Nut (see 21CFR, Sec. 888.3030)

4. **Device Classification for the subject device and/or predicate device:**

Class II

5. **Proposed Regulatory Class:** Class II

6. **Device Product Code:** HRS

7. **Subsequent Product Codes:** JDW, HWC, HTN

8. **Panel Code:** Orthopedics/Part 87

Section IV: Device Background

The I.T.S. Pelvic Reconstruction System (PRS) is a low profile - multiple design - pelvic plating, sacral rod securing and guided cannulated screw - group of fracture fixation systems - for fracture fixation of disruptive bone fracture injuries to the pelvic region in an adult or pediatric patient.

The I.T.S. PRS Low Profile – Multiple Type – Pelvic Plating System comes with a full complement of pelvic plate types in multiple hole and length sizes as: (1) Straight, (2) 108mm - radius Curved, (3) 88mm - radius Curved, (4) Symphysis, (5) SIJ L-shaped Left & Right, and (6) SIJ Closed. All plate types have a smooth and rounded profile and are contourable in three planes. All low profile plates accept a screw-to-plate head locking feature as a 5.9mm cancellous screw to create a stable, fixed angle, construct between the plate and screw(s). A standard 4.5mm cortical and 5.9mm cancellous compression screw is also available to facilitate reduction and create compression between the plate and bone. All screws are pre-drilling and self-tapping in design. The plate design is made from CP Titanium, a known biocompatible material, with a TIODIZE II surface treatment – which improves the material surface hardness and biocompatibility to adjacent tissues. The high strength 6-4 Alloyed Titanium cancellous and cortical screws allow a $\pm 20^\circ$ angle for positioning in relation to the plate so as to secure maximum purchase with bone. The low-profile plate design minimizes soft-tissue irritation in the pelvic region for the patient. A full compliment of instrumentation is available for use with the system.

510(k) Submission
I.T.S. Implantat-Technologie-Systeme GmbH
Pelvic Reconstruction System (PRS)

The I.T.S. PRS 7.3mm Cannulated SIJ Traction Screw & Washer is a complement of various length cannulated screws in a 7.3mm size with a variable thread length for guided fixation of pelvic bone fractures. All screws feature a cancellous thread form and are pre-drilling, self-tapping, and back-tapping (if removal is necessary) in design. All screws are used in conjunction with a 3.2mm guide wire and x-ray imaging for precise placement in bone across a pelvic bone fracture site. Washers (flat and curved) are available for thin cortex or osteopenic bone where the screw head may break through. The cannulated screw and washers are made from high strength 6-4 Alloyed Titanium material and have a TIODIZE II surface treatment for increased surface material hardness and biocompatibility to adjacent tissue. A full compliment of instrumentation is available for use with the system.

The I.T.S. PRS Sacral Rod System consists of a 'each-end threaded' pin, wedge-shaped washers, and nut/locknuts construct. One end of the threaded pin has a trocar point to assist in guiding the pin through pre-drilled holes if the guiding instrument is not used. The unit is assembled with one(1) 'each-end threaded' pin, two(2) wedge-shaped washers and two(2) nut/locknuts for stabilizing ilio-iliac posterior 'pelvic ring' disruption injuries. Typically, two assembled units are used for stabilizing the 'pelvic ring' injury. The threaded pin, wedge-shaped washers, and nut/locknut are manufactured from high strength 6-4 Alloyed Titanium material and are surface treated with the TIODIZE II surface treatment for increased material surface hardness and biocompatibility to adjacent tissues. Again, a full complement of instrumentation is available for use with the system.

510(k) Submission
I.T.S. Implantat-Technologie-Systeme GmbH
Pelvic Reconstruction System (PRS)

The attached brochure documents from I.T.S. GmbH describe the surgical technique and give examples of use for the multiple type - pelvic plating, sacral rod, and cannulated screw Pelvic Reconstruction System (PRS) along with a device listing of all implants and instrumentation.

EMERGENCY TEAM FOR BROKEN BONES


Pelvic Reconstruction System

PRS
 PELVIC RECONSTRUCTION SYSTEM



SACRAL RODS

In combination with I.T.S. sacral rods (7.3 mm diameter, 100 mm length, 1.5 mm thread) and I.T.S. spacer disks (flat & curved), there is now a complete solution for all established operative methods concerning injuries of the pelvic ring and acetabulum.



STRAIGHT PLATE

SACROILIAC - JOINT PLATE

SACRAL RODS

CURVED PLATE

SYMPHYSIS PLATE



A complete system for the treatment of all injuries to the pelvis and acetabulum

26

■ Pelvic Reconstruction System


Complete system for the treatment for all injuries of the pelvis and acetabulum

In combination with **I.T.S. sacral rods**, **I.T.S. 7.3 mm cannulated traction screws** [with adjustable thread] and **I.T.S. spacer disks** [flat / curved], there is now a complete solution for all established operative methods concerning injuries of the pelvic ring and acetabulum.

■ Indications

- All pelvic injuries

■ Combination - Complete Solution



I.T.S. IMPLANTAT-TECHNOLOGIE-SYSTEME GMBH IS RESPONSIBLE FOR THE CONTENTS AND MEDICAL FACTS.

Pelvic Reconstruction System

27

Pelvic Reconstruction System




SYMPHYSIS PLATE
PROFILE HEIGHT: 3.5 MM



SIJ-PLATE [CLOSED]
PROFILE HEIGHT: 2.5 MM
SIJ-PLATE [L-SHAPE]
PROFILE HEIGHT: 2.5 MM



CURVED PLATE [RADIUS 88 MM]
PROFILE HEIGHT: 2.5 MM




CURVED PLATE [RADIUS 108 MM]
PROFILE HEIGHT: 2.5 MM



STRAIGHT PLATE
[DORSAL ILIO - ILIAC DISTANCE OSTEOSYNTHESIS]
PROFILE HEIGHT: 3.5 MM

Characteristics



- Plate curved in 2 different degrees - anatomically contoured [for male and female pelvic anatomies]
- Can be anatomically shaped [bending pliers and setting instruments]
- Locking screw arrangement
- Bending templates [to avoid repeated in situ attachment and readjustment]
- Variable angle [$\pm 15^\circ$] between plate and screw - also in the case of locking screw arrangement
- Material of plate: pure titanium
- Material of screws: TiAL6V4 ELI
- I.T.S. surface treatment [will not grow into bone]
- 3 types of screws can be used



4.5 MM CORTICAL SCREW
5.9 MM CANCELLOUS-BONE SCREW
5.9 MM ANGLE-STABLE CANCELLOUS-BONE SCREW






Application of Acetabular Plate

- Access by means of Kocher-Langenbeck approach, possibly using trochanter-flip osteotomy and surgical luxation.
- Reduction and temporary fixation using K-wire.
- Shaping of the bending template.
- Appropriate final adjustment of the plate.
- Application of the plate and temporary fixation using spikes and fluoroscopic or X-ray inspection.
- Finally, insertion of cortical, cancellous-bone and locking screws into the holes of the plate.
- Fluoroscopic or X-ray inspection.
- Drainage of the area and layered closure of the wound [refixation of trochanter major if required].





■ Reconstruction of Anterior Pelvic Ring using Reconstruction Plate

- Ilio-inguinal approach or possibly Stoppa approach.
- Reduction using Weber or Jungbluth forceps, Schanz screw and temporary K-wire fixation if required.
- Attachment and adjustment of appropriate bending template, with fluoroscopic or X-ray inspection if required.
- Adjustment and setting of reconstruction plate on the bending template using setting instruments or bending pliers.
- Attachment of the plate and temporary fixation using spikes and fluoroscopic or X-ray inspection.
- Finally, insertion of cortical, cancellous-bone and locking screws. Conclude with fluoroscopic or X-ray inspection.
- Drainage, closure of wound.





Application of Symphysis Plate

- Lower medial laparotomy [emergency] or Pfannenstiel incision.
- Extra-peritoneal opening of the pelvis minor in the linea alba.
- Crenulation of attachments to the rectus abdominis from inside, if possible do not sever laterally.
- Resection using Weber or Jungbluth forceps.
- Positioning of 4 or 6-hole plate superiorly [adjust using setting instruments if required].
- Temporary fixation using spikes and fluoroscopic or X-ray inspection.
- Finally, insertion of cortical, cancellous-bone and locking screws.
- Drainage of cavum retzii, closure of the wound.







Application of SIJ Plate

- Antero-lateral approach or first window in the context of an ilio-inguinal approach.
- Adjustment of sacroiliac joint with Hohmann retractors.
- Resection using Weber or Jungbluth forceps.
- Adjustment of a double-hole plate or double L-plates.
- Temporary fixation using spikes and fluoroscopic or X-ray inspection.
- Finally, insertion of cortical, cancellous-bone and locking screws.
- Drainage, closure of wound.



Ilio-Iliac Distance Osteosynthesis



- Approach:
Bilateral superior 5-cm-long incision from the spina iliaca posterior superior.
- Bilateral exposure of rear iliac crest.
- Reduction using Schanz screw, longitudinal traction on leg, if necessary with the help of reduction forceps.
- Determination of plate length.
- Chisel off plate bearing [approx. 4 mm deep].
- Bend the plate from one side. Push the plate through behind the sacrum. Bend the plate on both sides in situ.
- Adjust plate and insert both screws on iliac crest. Alternate tightening of screws.
- Insert locking screws in both holes of the short side piece of the plate, at the same time making sure that the screwing angle is 10-15° to the plane of the already attached screw on iliac crest [to prevent collision of the screws - this is not a problem when using **I.T.S. locking screws**].
- Drainage, closure of wound.




33

SIJ Screw Fixation

- Approach:
Abdominal position, lateral stab incision after fluoroscopic inspection.
- Position the guide wire under fluoroscopic inspection, alternating in inlet and outlet view or with 2D or 3D navigation.
- Determine the length and position cannulated self-drilling 7.3 mm screw with spacer disk. [It is necessary to drill into cortical bone only when bone is very hard].
- Final X-ray inspection and wound closure.




Implants & Instruments




Instruments

	Item no.
A Bending forceps	66251-1
B Counter bending forceps	66251-2
C Plate benders	66252-14
D Sterilisation case	50168



35



Instruments

	Item no.
1 Ø 4.5 mm cortical screw 16-60 mm in 4 mm steps 60-80 mm in 5 mm steps	32455-xx
2 Ø 5.9 mm cancellous screw 16-60 mm in 4 mm steps 60-80 mm in 5 mm steps	30591-xx
3 Ø 5.9 mm cancellous stabilisation screw 16-60 mm in 4 mm steps 60-80 mm in 5 mm steps	37592-xx
4 Handle with AO - coupler	53011
5 Socket spanner SW 7 with handle	560701-350
6 Length gauge	59022
7 Ø 3.5mm drill L=280 mm	61353-280
8 Ø 3.5mm drill L=28mm	61353-110
9 Drill guide 2.5 / 3.5 mm	62252
10 Screwdriver shank 3.5 mm selfwedging AO	54353-180SH
11 Screwdriver shank 3.5 mm selfwedging AO	54353-90SH
12 Fixation screws wrench size 7	70301-7
13 Sterilisation case	50169


Implants

	Item no.
I Straight plate	21181-xx
II Curved plate R 108 mm	21195-xx
III Curved plate R 88 mm	21194-xx
IV Symphysis plate	21161-xx
V Sacroiliac-joint-plate right	21171-xx
Sacroiliac-joint-plate left	21172-xx
VI Sacroiliac-joint-plate closed	21173-xx
VII Sterilisation case	50167

Bending templates

	Item no.
Bending templates R 88 mm	67194-xx
Bending templates R 108 mm	67195-xx

36



IMPLANT-TECHNOLOGY-SYSTEMS

Item no.
32455-x
32591-x
32592-x
32611
32620-1-350
32622
32623-350
32624-110

- 1 6-4.5 mm cortical screw
- 1 6-6 mm in 4 mm steps 60-80 mm in 5 mm steps
- 2 6-5.9 mm cancellous screw
- 1 6-8 mm in 4 mm steps 60-80 mm in 5 mm steps
- 2 6-5.9 mm cancellous stabilization screw
- 1 6-6 mm in 4 mm steps 60-80 mm in 5 mm steps
- 4 Handle with AC-coupler
- 4 Socket spanner SW 7 with handle
- 2 Light gauge
- 1 6-3 mm drill 2-580 mm
- 1 6-3 mm drill 2-58 mm



I.T.S. GmbH
 AUTAL 28, A-8301 LASSNITZHÖHE / GRAZ
 AUSTRIA
 TEL +43 [0]316 211 21 0
 FAX +43 [0]316 211 21 20
 OFFICE@ITS-IMPLANT.COM
 WWW.ITS-IMPLANT.COM

© I.T.S. Spectromed Graz/Austria 2006. Errors, printing errors and technical changes reserved.

P.O. No. IPRS-0806-E CE 0297

I.T.S. Implante Technologie Systeme GmbH

Pelvic Reconstruction System (PRS) 7.3mm Cannulated Screw & Washer System

Surgical Technique

1. Prepare the patient with a general or regional anesthetic to the affected pelvic region area.
2. During the procedure, observe (using x-ray fluoroscopy) the fractured pelvic bone segment area(s). X-ray fluoroscopy observation is required during the entire procedure to ensure proper guide pin placement, proper bone screw depth, and adequate fracture impaction/reduction by the bone screw(s) device.
3. Following fracture reduction under image fluoroscopy control, insert the required size (3.2mm) calibrated guide pin across the pelvis fracture site engaging subchondral bone. If additional fixation is necessary across the fracture site, insert additional calibrated guide pins. Place the multiple guide pins parallel in direction to each other and in such a manner to avoid interference contact if washers are used.
4. Place the calibrated Depth Gauge over the calibrated guide pin and read the actual depth of the pin in the bone. The surgeon may elect to use the appropriate screw 5-10mm less than the Depth Gauge reading depending on screw penetration (ie.; so as not to penetrate into soft tissue).
5. If the cortex is thin or damaged and/or osteoporotic, a flat or curved washer may be placed under the head of the screw to prevent the screw head from being pulled into the cortex. A washer will distribute stress over the cortex of the bone.
6. Using the Gauge 5.0mm Cannulated Screwdriver, drive and seat the self-drilling and self-tapping screw over the guide pin and check the screw depth and fracture impaction/reduction with fluoroscopic x-ray.
7. Remove the calibrated guide pin.
8. Repeat steps 4 -7 for additional screw fixation across additional pelvic fracture sites and close.

The operative management of pelvic ring fractures has always represented a challenge to the treating surgeon. The sufficient stabilization of such an injury requires both a profound knowledge of the anatomy and biomechanics of the pelvis and an understanding of the mechanism of fracture. The chosen technique of a posterior stabilization depends on the one hand on the kind of fracture, and on the other hand on the personal experience of the individual operator.

Threaded rods as fixators of posterior pelvic ring instability have been well known and well proven for a long time, especially for sacral fractures in all three zones and for detachments of the SI joint. The technique was first described in 1913 by A. Lambotte, but frequent use of sacral rods only started in the 1980s. The stability of sacral rods has been proved in studies, and in many investigations it has been used as a "gold standard".

SACRAL-RODS

With the present technically perfected system, the technique has been standardized and considerably simplified. The use of titanium plates imposes no restrictions on postoperative imaging diagnosis.



For fixation of posterior pelvic ring instabilities



The operative management of **pelvic ring disruptions** has always represented a challenge to the treating surgeon. The sufficient stabilisation of such an injury requires both a profound knowledge of the anatomy and biomechanics of the pelvis and an understanding of the mechanism of fracture. The chosen technique of a posterior stabilisation depends on the one hand on the kind of fracture, and on the other hand on the personal experience of the individual operator.

Threaded rods as fixators of posterior pelvic ring instability have been well known and well proven for a long time, especially for sacral fractures in all three zones and for detachments of the SI joint (the technique was first described in 1913 by A. Lambotte, but frequent use of sacral rods only started in the 1980s.) The stability of sacral rods has been proved in numerous studies, and in many investigations the rods have been used as a „gold standard“ and reference value.

In the case of type C injuries, surgical intervention has first to be ventrally undertaken as a rule by means of plate osteosynthesis. The rear stabilisation is carried out either in the same session or in a quickly following second surgical intervention. The sacral rods are attached by means of the dorsal parts of the crista iliaca dorsal to the sacrum and dorsal to the sacro-iliac joint. The rods are secured on both sides by a nut and locknut.

With the present technically perfected armamentarium, the technique has been standardised and considerably simplified. The use of titanium plates imposes no restrictions on postoperative imaging diagnostics.




PROPERTIES

- Threaded rods with washers and locknuts for ilio-iliac compression osteosynthesis
- Stepped thread to increase fatigue strength
- Washer: wedge-shaped, ribbed, anatomically contoured
- Stable guiding instrument for simple and safe positioning of threaded rods
- Socket spanner system with AO standard adapter for fast and simple fastening of nuts
- Material: TiAl6V4 ELI
- I.T.S. surface treatment



ADVANTAGES

- Simple and safe operative technique
- Low possibility of iatrogenic injury to neural structures
- Low intraoperative X-ray exposure




INDICATIONS

- Dorsal stabilisation on the posterior pelvic ring for unstable pelvic ring injuries type C/AO,

especially in the case of:

- Unstable sacral fractures in all three zones
- Osteo-ligament instabilities of the sacroiliac joint



CONTRAINDICATIONS

- Fractures of the ala of the ilium in the posterior area
- Fractures of the os sacrum in the transforaminal area with CT-verified fractured parts and fragments requiring an open revision and decompression
- For the case that the relevant anatomy appears to preclude a safe positioning of the rods behind the sacral canal (rare)

DESCRIPTION OF SURGERY

VENTRAL STABILISATION

As a rule, a ventral stabilisation is carried out as the first intervention since a secondary reposition of the anterior pelvic ring in the case of an already reduced dorsal fracture could lead to an uncontrolled compression of the plexus sacralis.




DESCRIPTION OF SURGERY

DORSAL STABILISATION

Dorsal stabilisation is carried out in the prone position of the patient, usually under full anaesthesia and in the same session.


A curved skin incision is made on both sides lateral to the sacro-iliac joint, from the spina iliaca posterior superior extending cranially with a length of about 5 cm.



PREPARATION

Preparation of the subcutaneous fat tissue and carried down to the fascia of the glutaeus maximus muscle.

Loosen the musculature from the outer surface of the dorsal os ileum.




DESCRIPTION OF SURGERY

REPOSITION



A possibly necessary reposition of half of the pelvis is carried out by means of a Schanz screw attached to the crista iliaca.

Deployment of the guiding instrument and checking by means of image converter.



Guiding instrument:

It is strongly advised **not** to use the enclosed guiding device for reducing bone fragments since it could be damaged in this way and prevent a trouble-free attachment of the sacral rods.





DESCRIPTION OF SURGERY

The **linea glutea posterior** acts as the reference point.


The guiding instrument is positioned on a considered line between the origin of the linea glutea posterior on the iliac crest and the spina iliaca inferior; the drill hole will then lie some 0.5 cm higher.

The **caudal screw** is attached first, and the point of entry is to be chosen such that the rod is at a tangent to the crista centralis of the os sacrum or even perforates it, but at least comes to rest behind the lamina dorsalis of the sacrum.




The **guiding instrument** can be swivelled in the caudal direction during radioscopy (otherwise the guiding instrument can become loose due to divergence with the ala of the ileum).

The **threaded rod** is equipped with a drill bit, and young bones should be predrilled. After that, the threaded rod is introduced over the guiding instrument.



Pre-drilling:

For predrilling we recommend the use of the enclosed drill to achieve a satisfying results.




After the **threaded rod** has penetrated the second corticalis of the opposite os ileum, the guiding instrument is dismantled on this side.




45


DESCRIPTION OF SURGERY



The asymmetric **washer, nut and locknut** are attached (this is done the most simply by continuing to turn the rod with the drill). Both nuts must be definitively tightened.



The drill chuck is released, the guiding instrument removed, and the washer is attached as are, successively, both nuts on the second side using the special socket spanner and the soft-part protector.





Wedge-shaped spacers:

Please note the exact anatomical position shown here to arrive at the best result.

The threaded rod is equipped with a drill bit and young bones should be drilled. After that the threaded rod is introduced over the guiding instrument.

Pre-drilling: For pre-drilling we recommend the use of the enclosed drill to achieve a satisfying result.





46

DESCRIPTION OF SURGERY

The first nut is tightened until the desired compression of the fracture area is reached (care must be taken not to use excessive compression with transforaminal fractures of the sacrum.) Both nuts must be tightened. Remove the excess end using the bolt cutters.

The same procedure is undertaken with the second screw, which is to be positioned some 2 cm cranially to the first screw.

**CHECK**

Final check using image converter or X-ray, redon drainage, subsequent closure of the wound.

POSTOPERATIVE TREATMENT

Dependent on type of fracture:


In the case of instabilities of the SI joint and transalar fractures, we recommend a partial stress of the injured side for six weeks, then moving on to full stress.

In the case of transforaminal and central fractures, we recommend the injured side not be subject to stress for six weeks, then six weeks partial stress, moving on to full stress after CT check after the 12th week at the earliest.

The materials should be removed since two joints are bridged over.

Removal of the implants is carried out as a rule after one year at the earliest.

IMPLANTS & INSTRUMENTS



48




...

	item no.
1 Threaded Rod	39401-01
2 Washer, wedge-shaped, ribbed (in round sieve)	39401-02
3 Counternut , gauge 10 (in round sieve)	39401-03
4 Counternut , gauge 7 (in round sieve)	39401-04

INSTRUMENTS

	item no.
5 Guiding Instrument	62403
6 Socket Spanner, gauge 7, l=240	56702-240
7 Socket Spanner, gauge 10, l=200	561002-200
8 Tissue Protection Sleeve	62160-180
9 Spiral Drill d=4mm, l=225mm	61403-225
10 Flat spanner, gauge 7	70007
11 Flat spanner, gauge 10	70010
12 AO - Screw Adapter	53400
13 Bolt Cutter	65394
14 Sterilisation Case	50140



49



Item no.
39401-01
39401-02
39401-03
39401-04

4. Counter nut, gauge 7 (in round sieve)
3. Counter nut, gauge 10 (in round sieve)
2. Washer, wedge-shaped, tipped (in round sieve)
1. Threaded Rod



Item no.
62403
66703-240
661003-500
62160-160

8. Tissue Protection Sleeve
7. Socket Spanner, gauge 10, F=200
6. Socket Spanner, gauge 7, F=240
5. Guiding instrument

I.T.S. GmbH
AUTAL 28, A-8301 LASSNITZHOHE / GRAZ
TEL +43 (0)316 211 21 0
FAX +43 (0)316 211 21 20
OFFICE@ITS-IMPLANT.COM
WWW.ITS-IMPLANT.COM

© I.T.S. Spectromed Graz/Austria 2005. Errors, printing errors and technical changes reserved.

P.O.No. I-SS-0305-E CE 0297

50

I.T.S. Implantat Technologie Systeme GmbH

Pelvic Reconstruction System (PRS)

Implant and Instrument Listing

PRS Sacral Rod:

IMPLANTS

<u>Part Number</u>	<u>Description</u>
39401-01	Threaded Rod
39401-02	Washer, wedge-shaped, ribbed
39401-03	Counter-nut, gauge 10
39401-04	Counter-nut, gauge 7

INSTRUMENTS

62403	Guiding Instrument
56702-240	Socket Spanner, gauge 7, L=240
561002-200	Socket Spanner, gauge 10, L=200
62160-180	Tissue Protection Sleeve
61403-225	Spiral Drill d=4mm, L=225mm
70007	Flat Spanner, gauge 7
70010	Flat Spanner, gauge 10
53400	AO – Screw Adapter
65394	Bolt Cutter
50140	Sterilization Case

PRS 7.3mm Cannulated SIJ Traction Screw: - Previously Approved under K060156

IMPLANTS

<u>Part Number</u>	<u>Size</u>	<u>Description</u>
31734-35	7.3mm x 35mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 18mm Variable Th'd Length
31734-40	7.3mm x 40mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 20mm Variable Th'd Length
31734-45	7.3mm x 45mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 23mm Variable Th'd Length
31734-50	7.3mm x 50mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 25mm Variable Th'd Length
31734-55	7.3mm x 55mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 28mm Variable Th'd Length
31734-60	7.3mm x 60mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 30mm Variable Th'd Length
31734-65	7.3mm x 65mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 33mm Variable Th'd Length

Implant and Instrument Listing Continued:

<u>Part Number</u>	<u>Size</u>	<u>Description</u>
31734-70	7.3mm x 70mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 35mm Variable Th'd Length
31734-75	7.3mm x 75mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 38mm Variable Th'd Length
31734-80	7.3mm x 80mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 40mm Variable Th'd Length
31734-85	7.3mm x 85mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 43mm Variable Th'd Length
31734-90	7.3mm x 90mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 45mm Variable Th'd Length
31734-95	7.3mm x 95mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 48mm Variable Th'd Length
31734-100	7.3mm x 100mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 40mm Variable Th'd Length
31734-105	7.3mm x 105mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 40mm Variable Th'd Length
31734-110	7.3mm x 110mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 40mm Variable Th'd Length
31734-115	7.3mm x 115mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 40mm Variable Th'd Length
31734-120	7.3mm x 120mm	7.3mm Cannulated Cancellous Screw, 7.3mmØ, 40mm Variable Th'd Length

PRS Washer

36731	7.5mm x 15mm	Washer, Flat, for 7.3mm Cannulated Cancellous, 7.5mm ID x 15.0mm OD
36732	7.5mm x 14mm	Washer, Curved, for 7.3mm Cannulated Cancellous, 7.5mm ID x 14.0mm OD

INSTRUMENTS

59321	Depth Gauge for calibrated Guide Wire
35324-228	3.2mm diameter Guide Wire, Stainless Steel, Threaded Tip, L=228mm
54502-120	Screwdriver, cannulated, Allen Socket Insert, 5.0mm gauge, L=120mm
53032	Handle for Screwdriver, Flat sides
50178	Sterilization Case

PRS Plating and Screw Fixation Systems:

PRS Straight Plate IMPLANT

<u>Part Number</u>	<u>Size</u>	<u>Description</u>
21181-10	10	Straight Plate, 10 Hole
21181-11	11	Straight Plate, 11 Hole
21181-12	12	Straight Plate, 12 Hole
21181-13	13	Straight Plate, 13 Hole
21181-14	14	Straight Plate, 14 Hole

PRS Curved Plate, R=108mm, IMPLANT

<u>Part Number</u>	<u>Size</u>	<u>Description</u>
21195-4	4/R108	Curved Plate, 4 Hole, Radius = 108mm
21195-5	5/R108	Curved Plate, 5 Hole, Radius = 108mm
21195-6	6/R108	Curved Plate, 6 Hole, Radius = 108mm
21195-7	7/R108	Curved Plate, 7 Hole, Radius = 108mm
21195-8	8/R108	Curved Plate, 8 Hole, Radius = 108mm
21195-10	10/R108	Curved Plate, 10 Hole, Radius = 108mm
21195-12	12/R108	Curved Plate, 12 Hole, Radius = 108mm
21195-14	14/R108	Curved Plate, 14 Hole, Radius = 108mm
21195-16	16/R108	Curved Plate, 16 Hole, Radius = 108mm

PRS Curved Plate, R=88mm, IMPLANT

21194-4	4/R88	Curved Plate, 4 Hole, Radius = 88mm
21194-5	5/R88	Curved Plate, 5 Hole, Radius = 88mm
21194-6	6/R88	Curved Plate, 6 Hole, Radius = 88mm
21194-7	7/R88	Curved Plate, 7 Hole, Radius = 88mm
21194-8	8/R88	Curved Plate, 8 Hole, Radius = 88mm
21194-10	10/R88	Curved Plate, 10 Hole, Radius = 88mm
21194-12	12/R88	Curved Plate, 12 Hole, Radius = 88mm
21194-14	14/R88	Curved Plate, 14 Hole, Radius = 88mm
21194-16	16/R88	Curved Plate, 16 Hole, Radius = 88mm

PRS Symphysis Plate IMPLANT

21161-4	4	Symphysis Plate, 4 Hole
21161-6	6	Symphysis Plate, 6 Hole

PRS SacroIliac-Joint L-Shape IMPLANT

21171-5	Left	SIJ-plate, L-shape, 5 Hole, Left
21172-5	Right	SIJ-plate, L-shape, 5 Hole, Right

PRS SacroIliac-Joint Closed IMPLANT

21173-4	Closed	SIJ-plate, Closed, 4 Hole
---------	--------	---------------------------

Implant and Instrument Listing Continued:

PRS 5.9mm Cancellous Bone Screw, Angle Stable, IMPLANT

<u>Part Number</u>	<u>Size</u>	<u>Description</u>
37592-16	5.9mm x 16mm	5.9mm Cancellous Bone Screw, Angle Stable, 16mm Length
37592-20	5.9mm x 20mm	5.9mm Cancellous Bone Screw, Angle Stable, 20mm Length
37592-24	5.9mm x 24mm	5.9mm Cancellous Bone Screw, Angle Stable, 24mm Length
37592-28	5.9mm x 28mm	5.9mm Cancellous Bone Screw, Angle Stable, 28mm Length
37592-32	5.9mm x 32mm	5.9mm Cancellous Bone Screw, Angle Stable, 32mm Length
37592-36	5.9mm x 36mm	5.9mm Cancellous Bone Screw, Angle Stable, 36mm Length
37592-40	5.9mm x 40mm	5.9mm Cancellous Bone Screw, Angle Stable, 40mm Length
37592-44	5.9mm x 44mm	5.9mm Cancellous Bone Screw, Angle Stable, 44mm Length
37592-48	5.9mm x 48mm	5.9mm Cancellous Bone Screw, Angle Stable, 48mm Length
37592-52	5.9mm x 52mm	5.9mm Cancellous Bone Screw, Angle Stable, 52mm Length
37592-56	5.9mm x 56mm	5.9mm Cancellous Bone Screw, Angle Stable, 56mm Length
37592-60	5.9mm x 60mm	5.9mm Cancellous Bone Screw, Angle Stable, 60mm Length
37592-65	5.9mm x 65mm	5.9mm Cancellous Bone Screw, Angle Stable, 65mm Length
37592-70	5.9mm x 70mm	5.9mm Cancellous Bone Screw, Angle Stable, 70mm Length
37592-75	5.9mm x 75mm	5.9mm Cancellous Bone Screw, Angle Stable, 75mm Length
37592-80	5.9mm x 80mm	5.9mm Cancellous Bone Screw, Angle Stable, 80mm Length

PRS 5.9mm Cancellous Bone Screw IMPLANT

<u>Part Number</u>	<u>Size</u>	<u>Description</u>
30591-16	5.9mm x 16mm	5.9mm Cancellous Bone Screw, 16mm Length
30591-20	5.9mm x 20mm	5.9mm Cancellous Bone Screw, 20mm Length
30591-24	5.9mm x 24mm	5.9mm Cancellous Bone Screw, 24mm Length
30591-28	5.9mm x 28mm	5.9mm Cancellous Bone Screw, 28mm Length
30591-32	5.9mm x 32mm	5.9mm Cancellous Bone Screw, 32mm Length
30591-36	5.9mm x 36mm	5.9mm Cancellous Bone Screw, 36mm Length
30591-40	5.9mm x 40mm	5.9mm Cancellous Bone Screw, 40mm Length
30591-44	5.9mm x 44mm	5.9mm Cancellous Bone Screw, 44mm Length
30591-48	5.9mm x 48mm	5.9mm Cancellous Bone Screw, 48mm Length
30591-52	5.9mm x 52mm	5.9mm Cancellous Bone Screw, 52mm Length
30591-56	5.9mm x 56mm	5.9mm Cancellous Bone Screw, 56mm Length
30591-60	5.9mm x 60mm	5.9mm Cancellous Bone Screw, 60mm Length
30591-65	5.9mm x 65mm	5.9mm Cancellous Bone Screw, 65mm Length
30591-70	5.9mm x 70mm	5.9mm Cancellous Bone Screw, 70mm Length
30591-75	5.9mm x 75mm	5.9mm Cancellous Bone Screw, 75mm Length
30591-80	5.9mm x 80mm	5.9mm Cancellous Bone Screw, 80mm Length

89

Implant and Instrument Listing Continued:

PRS 4.5mm Cortical Bone Screw IMPLANT

<u>Part Number</u>	<u>Size</u>	<u>Description</u>
32455-16	4.5mm x 16mm	4.5mm Cortical Bone Screw, Std., 16mm Length
32455-20	4.5mm x 20mm	4.5mm Cortical Bone Screw, Std., 20mm Length
32455-24	4.5mm x 24mm	4.5mm Cortical Bone Screw, Std., 24mm Length
32455-28	4.5mm x 28mm	4.5mm Cortical Bone Screw, Std., 28mm Length
32455-32	4.5mm x 32mm	4.5mm Cortical Bone Screw, Std., 32mm Length
32455-36	4.5mm x 36mm	4.5mm Cortical Bone Screw, Std., 36mm Length
32455-40	4.5mm x 40mm	4.5mm Cortical Bone Screw, Std., 40mm Length
32455-44	4.5mm x 44mm	4.5mm Cortical Bone Screw, Std., 44mm Length
32455-48	4.5mm x 48mm	4.5mm Cortical Bone Screw, Std., 48mm Length
32455-52	4.5mm x 52mm	4.5mm Cortical Bone Screw, Std., 52mm Length
32455-56	4.5mm x 56mm	4.5mm Cortical Bone Screw, Std., 56mm Length
32455-60	4.5mm x 60mm	4.5mm Cortical Bone Screw, Std., 60mm Length
32455-65	4.5mm x 65mm	4.5mm Cortical Bone Screw, Std., 65mm Length
32455-70	4.5mm x 70mm	4.5mm Cortical Bone Screw, Std., 70mm Length
32455-75	4.5mm x 75mm	4.5mm Cortical Bone Screw, Std., 75mm Length
32455-80	4.5mm x 80mm	4.5mm Cortical Bone Screw, Std., 80mm Length

PRS Plating and Screw Fixation INSTRUMENTS

<u>Part Number</u>	<u>Size</u>	<u>Description</u>
67194-4	4/R88	Bending Template for Curved Plate, Radius = 88mm, 4 Hole
67194-5	5/R88	Bending Template for Curved Plate, Radius = 88mm, 5 Hole
67194-6	6/R88	Bending Template for Curved Plate, Radius = 88mm, 6 Hole
67194-7	7/R88	Bending Template for Curved Plate, Radius = 88mm, 7 Hole
67194-8	8/R88	Bending Template for Curved Plate, Radius = 88mm, 8 Hole
67194-10	10/R88	Bending Template for Curved Plate, Radius = 88mm, 10 Hole
67194-12	12/R88	Bending Template for Curved Plate, Radius = 88mm, 12 Hole
67194-14	14/R88	Bending Template for Curved Plate, Radius = 88mm, 14 Hole
67194-16	16/R88	Bending Template for Curved Plate, Radius = 88mm, 16 Hole
67195-4	4/R108	Bending Template for Curved Plate, Radius = 108mm, 4 Hole
67195-5	5/R108	Bending Template for Curved Plate, Radius = 108mm, 5 Hole
67195-6	6/R108	Bending Template for Curved Plate, Radius = 108mm, 6 Hole
67195-7	7/R108	Bending Template for Curved Plate, Radius = 108mm, 7 Hole
67195-8	8/R108	Bending Template for Curved Plate, Radius = 108mm, 8 Hole
67195-10	10/R108	Bending Template for Curved Plate, Radius = 108mm, 10 Hole
67195-12	12/R108	Bending Template for Curved Plate, Radius = 108mm, 12 Hole
67195-14	14/R108	Bending Template for Curved Plate, Radius = 108mm, 14 Hole
67195-16	16/R108	Bending Template for Curved Plate, Radius = 108mm, 16 Hole
66251-1		Bending Forceps
66251-2		Counter Bending Forceps
66252-14		Plate Benders
53011		Handle
560701-350		Socket Spanner SW 7 with Handle
59022		Length Gauge

55

PRS Plating and Screw Fixation INSTRUMENTS Continued:

<u>Part Number</u>	<u>Size</u>	<u>Description</u>
61353-280		3.5mm Drill, L = 280mm
61353-110		3.5mm Drill, L = 110mm
62252		Drill Guide 2.5/3.5mm
54353-180SH		Screwdriver SW 3.5mm
54353-90SH		Screwdriver SW 3.5mm
70301-7		Fixation Screws SW 7
50167		Sterilization Case (for Implants)
50168		Sterilization Case (for Bending Forceps, etc.)
50169		Sterilization Case (for Instrumentation)

510(k) Submission
I.T.S. Implantat-Technologie-Systeme GmbH
Pelvic Reconstruction System (PRS)

Section V: Intended Uses/Indications

The **intended use** of the I.T.S. Pelvic Reconstruction System (PRS) is to stabilize one or more pelvic bone fractures in the pelvic ring area of the pelvis of an adult or pediatric patient which include the use of the following pelvic fracture fixation systems.

Indications for use of the I.T.S. PRS Low Profile - Multiple Type - Pelvic Plating System include:

- 1) Fracture reconstruction of the acetabulum, sacrum, ilium, and entire pelvic ring,
- 2) Revision surgery of pseudoarthroses, non-unions and mal-unions,
- 3) Ilio-Iliac distance osteosynthesis,
- 4) Osteotomies,
- 5) Arthrodesis,
- 6) Sacroiliac joint dislocations, and
- 7) Symphysis pubis disruptions

Indications for use of the I.T.S. PRS Sacral Rod System include:

- 1) Fixation of fractures of the posterior pelvis,
- 2) Fixation of fractures of the posterior iliac spine,
- 3) Fixation of fractures of the posterior inferior iliac spine,
- 4) Dorsal stabilization of the posterior pelvic ring for unstable pelvic ring injuries,
- 5) Fixation of sacral fractures, and
- 6) Fracture dislocations of the sacro-iliac joint

Indications for use of the I.T.S. PRS 7.3mm Cannulated SIJ Traction Screw & Washer include fracture fixation of pelvic bone fractures where indicated and for periacetabular osteotomies.

for Δ see Section 13. RPJ 12/22/2006

~~The system is not intended for spinal use.~~

Contraindications for use include active infection near the fracture site, skeletally immature patients, severe osteoporosis or inadequate bone stock, fractures that cannot be acceptably reduced by traction alone, foreign body (material) sensitivity, long delay after injury, and patients with mental or neurological conditions who are unwilling or incapable of following postoperative care instructions.

510(k) Submission
I.T.S. Implantat-Technologie-Systeme GmbH
Pelvic Reconstruction System (PRS)

Contraindications for the sacral rod also include:

- 1) Fractures of the ala of the ilium in the posterior area,
- 2) Fractures of the os sacrum in the transforaminal area with CT-verified fractured parts and fragments requiring an open revision and decompression, and
- 3) For the case that the relevant anatomy appears to preclude a safe positioning of the rods behind the sacral canal.

Section VI: Device Description

The I.T.S. Pelvic Reconstruction System (PRS) encompasses a number of fracture fixation subsystems (multiple pelvic plate designs, sacral threaded rod, and cannulated screw & washer) to make available the necessary implant hardware to the physician in fracture fixation and reconstruction of pelvic ring fractures in the pelvis.

The I.T.S. PRS Low Profile – Multiple Type – Pelvic Plating System consists of the following plate types: 1) A Straight Plate at a 4.0mm thickness with 10 to 14 hole length sizes, 2) A Curved Plate with a 108mm radius at a 2.5mm thickness with 4 to 16 hole length sizes, 3) A Curved Plate with a 88mm radius at a 2.5mm thickness also with 4 to 16 hole length sizes, 4) A Symphysis Plate at a 4.0mm thickness in both a 4 & 6 hole size, 5) A Sacroiliac-Joint (SIJ) L-Shaped Plate at a 2.5mm thickness in a left and right 5 hole size, and 6) A SIJ Closed Plate at a 2.5mm thickness in a 4 hole size. All plate designs are low profile in thickness, made from CP Titanium for easy 3-Dimensional contouring to the pelvic anatomy, and are smooth and rounded to reduce any soft tissue irritation. The Pelvic Plating System also encompasses a number of cancellous (5.9mm Cancellous standard compression & 5.9mm Cancellous Locking Screws) and cortical (4.5mm standard compression screw) screw types and length sizes to accommodate the bone reconstruction and fixation requirements for various 'indications of use' with all multiple type Pelvic Plate designs. All bone screws are pre-drilling and self-tapping in design and manufactured from high strength 6-4 Alloyed Titanium. All components

510(k) Submission
I.T.S. Implantat-Technologie-Systeme GmbH
Pelvic Reconstruction System (PRS)

(plates & screws) have a TIODIZE II surface treatment to improve surface material hardness and enhance biocompatibility to adjacent tissue.

The design intent of PRS Pelvic Plating System is to introduce locking bone screws into both the bone and bone plate and provide stabilization to the fracture or osteotomy site. This concept prevents any loosening of the screw(s) from the plate and gives a solid construct of screw-to-plate across the fracture site. The bone screws are comprised of a harder and high strength 6-4 Alloy Titanium material where the outer head portion of the screw is configured with a 'tapered thread' and introduced (self-tapping) into the bone mass through the countersunk portion of the receiving hole in the softer material CP Titanium plate. The head-thread engages into the countersunk plate hole and locks into place within the plate at a pre-determined angle of up to 20°. This concept of the screw locking into the plate is comparable to other plating systems such as the I.T.S. Volar Radius Plate (K033756) and the Hand Innovations Volar Radius Plate (K002775).

The I.T.S. PRS Sacral Rod System consists of a threaded pin, wedge-shaped washer(2ea.), and nut/locknut(2ea.) design. One end of the threaded pin has a trocar point for guiding the pin if the guiding instrument is not used. The assembled unit is used to stabilize ilio-iliac posterior pelvic ring disruption injuries. All components are manufactured from high strength 6-4 Alloyed Titanium material and have the TIODIZE II surface preparation, again for increased material surface hardness and biocompatibility.

510(k) Submission
I.T.S. Implantat-Technologie-Systeme GmbH
Pelvic Reconstruction System (PRS)

The I.T.S. PRS 7.3mm Cannulated SIJ Traction Screw & Washer is a pre-drilling, self-tapping, and back-tapping (if removal is necessary) screw design for guided reduction of pelvic bone fractures. The system uses a 3.2mm calibrated guide wire for precise positioning of the cannulated screw across the pelvic bone fracture when used in conjunction with x-ray imaging. Washers (flat & curved) are available for use with the screw when encountering thin cortex or osteopenic bone – where the screw head may break through. All 7.3mm Cannulated Screws and Washers are manufactured from high strength 6-4 Alloyed titanium and have the TIODIZE II surface preparation for increased material surface hardness and biocompatibility.

All Pelvic Reconstruction fracture fixation Systems (PRS) have a full set of instrumentation for use with each system.

Marketing brochures are enclosed in this section to show this screw plate locking system design concept from the referenced 510k's. A surgical technique and listing of the implant system with instrumentation is enclosed in **Section IV** of this application.


Information regarding the I.T.S. PRS System is disclosed in the attached brochure (*entitled: PRS & Sacral - Rods*) within **Section IV** of this submission. The brochure describes in detail the Design Properties, Materials, Indications, and **listing** of *ALL size components for both implants and instrumentation*. Also, within Section IV is the surgical technique describing use of all devices.

Drawings of ALL implant components (with all sizes) within the system is enclosed in this section.



PROPERTIES


- Angular stability between screws and plate
- Free selection of screw angle to the normal of the plate
- Left and right versions
- Profile height: only 1.5 mm
- Anatomical shape
- Angle-stable 3.5mm spongiosa stabilization screw of 14 - 30mm, in 2mm steps
- Plate material: pure titanium
- Screw material: TiAl6V4
- Surface treatment




INDICATIONS

- Complex intra-articular fractures of the distal radius
- Complex extra-articular fractures of the distal radius
- Osteotomies of the distal radius


SURGICAL TECHNIQUE



VOLAR ACCESS via the tendon of the *flexor carpi radialis* muscle
Separation of the transverse carpal ligament.
Open reduction



FILLING OF THE BROKEN AREAS
with bone substitute (e.g. Tutobone®)



ADJUSTMENT OF THE PLATE
to the anatomy

Section XII: 510(k) Summary of Safety and Effectiveness

FEB - 6 2004

SAFE MEDICAL DEVICES ACT OF 1990
510(k) Summary

K033756
page 1 of 1

NAME OF FIRM: I.T.S. Implantat-Technologie-Systeme GmbH.
Autal 28.
Lassnitzhoche A - 8301
AUSTRIA

510(k) FIRM CONTACT: Al Lippincott
Engineering Consulting Services, Inc.
3150 E. 200th St.
Prior Lake, MN 55372

TRADE NAME: Volar Radius Plate with Angular Stability

COMMON NAME: Fracture Fixation Plating system for fracture fixation of the end of long bones

CLASSIFICATION: Single/multiple component metallic bone fixation appliances and accessories (see 21 CFR, Sec. 888.3030).

DEVICE PRODUCT CODE: HRS

SUBSTANTIALLY EQUIVALENT DEVICES: Avanta SCS/V Distal Radius Plate - Volar
Hand Innovations DVR Plate
Synthes Locking Distal Radius System 2.4

DEVICE DESCRIPTION: The I.T.S. Volar Radius Plate with Angular Stability provides various width 4 and 6 hole standard plates, various width 8, 10, and 12 hole long plates, various length stabilization screws, and various length cortical fixation screws. The volar radius plates are made from CP Titanium according to ASTM F 67-00 and the screws are made from 6-4 Alloyed Titanium according to ASTM F 136-98.

INTENDED USE: The I.T.S. Volar Radius Plate is used to stabilize distal radius fractures in the wrist with an accurate retention of articular fracture elements.


BASIS OF SUBSTANTIAL EQUIVALENCY: The I.T.S. Volar Radius Plate is substantially equivalent to the Avanta, Hand Innovations, and Synthes fracture plating systems

SUMMARY OF SAFETY AND EFFECTIVENESS: The I.T.S. Volar Radius Plate system is shown to be safe and effective for use in fracture fixation of the distal radius in the wrist.

THE ONLY VOLAR PLATE DESIGNED FOR DORSAL FRACTURES

PROXIMAL PEG ROW TO SUPPORT DORSAL ASPECT OF SUBCHONDRAL PLATE

PROXIMAL K-WIRE HOLES FOR TEMPORARY FIXATION TO PROXIMAL FRAGMENT




DISTAL PEG ROW TO SUPPORT VOLAR ASPECT OF SUBCHONDRAL PLATE

DISTAL FIXATION OPTIONS:

- SMOOTH PEGS OFFER THE STRONGEST SUPPORT
- THREADED PEGS TO LAG DORSAL FRAGMENTS
- CANCELLOUS SCREWS FOR VOLAR FRACTURES

DISTAL K-WIRE HOLES FOR TEMPORARY FIXATION AND PLATE ALIGNMENT TO DISTAL FRAGMENT



FUNCTIONAL RESULTS



8905 SW 87TH AVENUE, SUITE 220, MIAMI, FLORIDA 33176
 PHONE (305) 412-8010 (800) 800-8188 FAX (305) 412-8060
 WWW.HANDINNOVATIONS.COM

64

**Summary of Safety and Effectiveness As Required by 807.92(c)
for the
Distal Volar Radius Fracture Repair System**

DEC - 5 2000

Submitted by

Hand Innovations, Inc.
8905 SW 87th Avenue
Suite 100
Miami, FL 33176-2227
Phone: (305) 661-3054
Fax: (305) 595-6620

Contact Person: Al Weisenborn

Device Trade Name: Distal Volar Radius Fracture Repair System
Common Name: Distal Volar Radius Locking Plate with pegs and screws
Classification Name: Single/multiple Component Metallic Bone Fixation Appliances and Accessories, per § 888.3030

Identification of a Legally Marketed Predicate Device

The Hand Innovations, Inc. Distal Volar Radius Fracture Repair System (DVR) is substantially equivalent to the Distal Radius Plate System that is manufactured by manufactured and distributed by Synthes (USA) pursuant to K982732.

Device Description

The DVR system consists of a stabilization plate, four bone screws, and four fixation pegs. The screws are used to fix the proximal segment of the plate to the diaphysis and the pegs the distal bone fragment(s). The required components and specialized instruments are supplied non-sterile in a container suitable for steam sterilization. All components and instruments may be purchased separately.

Intended Use

The Distal Volar Radius Fracture Repair System is intended for the volar fixation of fractures and osteotomies involving the distal radius.

Summary of Technological Characteristics

The DVR is substantially equivalent to equivalent to the Distal Volar Radius Fracture Repair System that is manufactured by manufactured and distributed by Synthes (USA) pursuant to K982732. This has been demonstrated through a 15 point comparison of technological characteristics.

Summary of Performance Data

A dimensional analysis of the DVR fracture repair system components met design requirement. Tensile strength performance characteristics of the DVR and predicate device were tested; the DVR was found to meet or exceed the measured performance characteristics of the predicate device.

Conclusion

The DVR has been demonstrated to be equivalent to the Distal Radius Plate System that is manufactured and distributed by Synthes (USA) pursuant to K982732, Inc. by bench testing of both devices and comparison of technological characteristics.

The tissue/bone contact material of the device have been carefully selected for its long history of biocompatibility. The material meets the requirements of a recognized consensus standard for implantable orthopedic material.

The DVR was designed utilizing design controls compliant with the Quality System Regulation. The DVR will be manufactured per specifications and good manufacturing practices that ensure the device is safe and effective for its intended use.

PRS Sacral Rod System Drawings

CONFIDENTIAL

PRS 7.3mm Cannulated SIJ
Traction Screw System
Drawings

CONFIDENTIAL

PRS Pelvic Plating System Drawings

CONFIDENTIAL

Section VII: Materials

(b)(4)Product Specs



**510(k) Submission
I.T.S. Implantat-Technologie-Systeme GmbH
Pelvic Reconstruction System (PRS)**

Section VIII: Testing Information

Section VII in this submission produces the mechanical/chemical testing of the material as stipulated in the ASTM Specifications F 136 & F 67.

Enclosed in this section is an Engineering Rationale disclosing information as to why mechanical torque and dynamic bend testing is not necessary for the I.T.S. Pelvic Reconstruction System (PRS) due to similarities to predicate devices. Enclosed with this rationale are marketing references, dimensional and material information of the predicate device(s), and 510(k) Summary Statements of the predicate devices.

AUG 24 2000

3.0 Summary of Safety and Effectiveness Information

SPONSOR: Synthes (USA)
1690 Russell Road
Paoli, PA 19301
(610) 647-9700
Contact: Angela Silvestri

DEVICE NAME: Synthes Sacral Bar System

CLASSIFICATION: 21 CFR 888.3040
Smooth or threaded metallic bone fixation fastener.

PREDICATE DEVICES: Zimmer Threaded Sacral Rod
Synthes Threaded Bolt

DEVICE DESCRIPTION: The Synthes Sacral Bar System consists of a threaded bar, washers, and nuts. The bars are fully threaded. One end of the bar has a trocar point to guide the bar through pre-drilled holes. The bars are available in lengths ranging from 120 to 260 mm, in 10 mm increments. The washers that are used with this system are oval shaped and are designed to slide freely along the bars. Both rounded and straight nuts are provided with this system; the rounded nuts mate with the washers to create compression, while the straight nuts are then added to wedge against the rounded nuts to maintain compression.

INTENDED USE: The Synthes Sacral Bar System is intended for fixation of fractures of the posterior pelvis, in areas of the posterior superior iliac spine and posterior inferior iliac spine, for sacral fractures and fracture-dislocations of the sacro-iliac joint.


MATERIAL: Implant quality stainless steel.

121

129

Sacral Bar Kit

TECHNIQUE GUIDE





Sacral Bar Kit

Indications

Synthes Sacral Bars are indicated for fixation of the posterior pelvis, in areas of the posterior superior iliac spine, posterior inferior iliac spine, for sacral fractures, and fracture dislocations of the sacroiliac joint.

Features

- Sacral washer prevents rounded sacral nut from pressing into cortical bone
- Trocar tip for easy passage through soft tissues
- Fully threaded, in lengths from 160 mm–260 mm, in 20 mm increments
- Kits available sterile and nonsterile



Technique

This procedure can be performed manually using a Universal Chuck with T-Handle, or with power equipment.

Instruments and implants required:

- 2 Sacral Bar Kits [199.93x]
- 6.0 mm Drill Bit, Jacobs chuck, 195 mm [310.60J] (included in kit)
- Rod Cutter [388.72]
- Combination Wrench, 11 mm width across flats [321.16] *or*
- Ratchet Wrench, 11 mm width across flats [321.20]

Note: *In the case of a bilateral fracture, it is necessary to have at least one side of stable reduction using an iliosacral screw or posterior plate in addition to inserting sacral bars.*


The correct length of the sacral bar should extend 25–30 mm on either side of the posterior iliac crests (bar will be cut to size intraoperatively).

1 Make incision

Make slightly curved bilateral incisions just lateral to the posterior superior sacroiliac joint. The incision on the fracture side should be large enough to allow for direct visualization and reduction.

2 Reduce and stabilize fracture

Reduce the fracture using pointed reduction clamps from the spinous process of the sacrum to the outer aspect of the posterior tubercle, or from the opposite side using pelvic reduction forceps. After reduction, two Kirschner wires should be placed across the sacroiliac joint for provisional stabilization.



151

124 1

Technique (continued)

3 Drill




Drill a set of gliding holes using a 6.0 mm Drill Bit. The holes should be located so that the first bar is placed at the level of L5/S1 interspace as seen on C-arm or by palpation. The second set of holes should be at least 1.5 cm distal to the first. When drilling the gliding hole, make sure that the drill exits the inner aspect of the posterior tubercle above the sacral lamina. Care must be taken to ensure that there is enough bone in the posterior tubercle to hold the second bar. If not, use of an iliosacral screw or posterior plate will be necessary.

4 Insert 6.0 mm Threaded Sacral Bar

Drive the sharp trocar tip through the predrilled hole in the contralateral iliac spine. The sacral lamina must be seen, to ensure that the bar passes dorsal to it and does not enter the sacral canal; it may pass through the sacral spinous processes.


5 Apply sacral washers and sacral nuts

Place a 6.0 mm Sacral Washer [296.758] on each end of the sacral bar to prevent the sacral nuts from pressing into the bone. Mate the 6.0 mm Sacral Nuts, rounded [296.757] to the sacral washers. Tighten with an 11 mm wrench to obtain compression.



6 Apply hexagonal nuts

Add 6.0 mm Hexagonal Nuts [296.759] to maintain compression. With a hexagonal wrench, wedge the hexagonal nuts against the rounded sacral nuts.





7 Place second bar

Repeat steps 4–6 for insertion of the second sacral bar. Two bars must be used to avoid rotation.





8 Remove Kirschner wires




Technique (continued)

9 Trim Sacral Bars

Use the Rod Cutter [388.72] to cut the ends of the sacral bars flush with the hexagonal nuts.



Rod Cutter
[388.72]



Final construct

Product Information

6.0 mm Threaded Sacral Bar Kits[◇]

199.931	160 mm
199.932	180 mm
199.933	200 mm
199.934	220 mm
199.935	240 mm
199.936	260 mm

Note:

Construct described in this Technique Guide requires two Sacral Bar Kits.

Each kit contains:

296.76x	6.0 mm Threaded Sacral Bar, 1 ea.
296.757	6.0 mm Sacral Nuts, rounded, 2 ea.*
296.758	6.0 mm Sacral Washers, 2 ea.*
296.759	6.0 mm Hexagonal Nuts, 2 ea.*
310.60J	6.0 mm Drill Bit, Jacobs chuck, 195 mm*

6.0 mm Threaded Sacral Bars*

296.761	160 mm
296.762	180 mm
296.763	200 mm
296.764	220 mm
296.765	240 mm
296.766	260 mm

Additionally Available

150.16	Compact Air Drive II Set or
105.957	Power Drive Set
321.16	Combination Wrench, 11 mm width across flats
321.20	Ratchet Wrench, 11 mm width across flats
388.72	Rod Cutter (<i>necessary to complete procedure</i>)
393.10	Universal Chuck with T-Handle
398.86	Pelvic Reduction Forceps, with pointed-ball tips, 400 mm length, extra large, long jaws, speed lock
511.20	Oscillating Drill Attachment

[◇] Kits available nonsterile or sterile-packed. Add "S" to catalog number to order sterile Sacral Bar Kit.

* Individual components may be purchased separately, nonsterile only.

MAR 20 2006

P.1/2

K060156

Section XII: 510(k) Summary of Safety and Effectiveness

**SAFE MEDICAL DEVICES ACT OF 1990
510(k) Summary**

NAME OF FIRM: I.T.S. Implantat-Technologie-Systeme GmbH
Aural 28
Lassnitzhoehe A - 8301
AUSTRIA

510(k) FIRM CONTACT: Al Lippincott
Engineering Consulting Services, Inc.
3150 E. 200th St
Prior Lake, MN 55372

TRADE NAME: Straight Plate with Angular Stability & Screw System

COMMON NAME: Bone Plate & Screw System

CLASSIFICATION: Plate, Fixation, Bone (see 21 CFR, Sec. 888.3030),
Screw, Fixation, Bone (see 21 CFR, Sec. 888.3040),
Washer, Bolt, Nut (see 21 CFR, Sec. 888.3030)

DEVICE PRODUCT CODE: HRS

SUBSEQUENT PRODUCT CODE: HWC, HTN

SUBSTANTIALLY I.T.S. GmbH Clavicularplate with Angular Stability (K050852)
Smith & Nephew Peri-Loc Locking Bone Plates and Locking Bone
Screw System (K051735)

EQUIVALENT DEVICES Zimmer Periaricular Locking Plates and Screws (K051098)
Synthes 4.5mm LCP Straight Reconstruction Plates (K051986)
Acumed Congruent Plate System (K012655)
I.T.S. GmbH FR.O.H. Calcaneus Repair System (K051642)
Synthes Sterile 3.5mm and 4.0mm Cannulated Screw (K963192)
Zimmer/Pioneer Cannulated Screw System (K603496)
Synthes (USA) Spherical Washers (K052483)

DEVICE DESCRIPTION: The I.T.S. Straight Plate with Angular Stability is a low-profile 4, 6, or 8 hole plate with various length cortical and/or cancellous self-tapping stabilization locking and/or compression screws. The Straight Plate is made from CP titanium according to ASTM F 67-00 and all screws are made from 6-4 Alloyed Titanium according to ASTM F 136-02. The plate and screws are surface conditioned with a TIODIZE, Type II preparation.

The I.T.S. Screw System is a group of cannulated fracture fixation screws in various diameters of 4.0mm, 6.5mm, and 7.3mm and lengths. A complement of flat and spherical Washers are available with the system. All screws and washers are made from 6-4 Alloyed Titanium according to ASTM F 136-02 and are surface conditioned with a TIODIZE, Type II preparation.

129

K 060156

P 2/2

INTENDED USE:

The intended use of the I.T.S. Straight Plate with Angular Stability is to stabilize an osteotomy or fracture of small bones, long bones, the pelvis and the calcaneus in an adult or pediatric patient.

Indications for use include comminuted fractures, supercondylar fractures, intra-articular and extra-articular condylar fractures, fractures in osteopenic bone, nonunions, and malunions. And, as well, a fracture or osteotomy of the tibia, fibula, femoral condyle, acetabulum, humerus, ulna, middle hand and middle foot bones; treatment of the calcaneal; hip arthrodesis, and provisional hole fixation. The 3.5mm Cortical and Cancellous Angle Stable Screws used in conjunction with the Straight Plate may be used only on small bones.

The intended use of the I.T.S. Screw System is for corrective osteotomy or internal fracture fixation of the patella, pelvis, ankle, and long bones in an adult or pediatric patient.

For the 4.0mm Cannulated Cancellous Screw, indications for use are for radial and ulnar fractures, fractures of the proximal/distal humerus and of the patella, and for tendon fixation, maisonneuve injuries and disruption of the syndesmosis with bimalleolar or supramalleolar fractures and the instability of the talus centering.

For the 6.5mm Cannulated Cancellous Screw, indications for use are for fractures of the femoral neck, tibiaplateau, of the sacrum and the articular cavity of the hip joint, and the metaphyseal fractures of the distal femur and distal tibia, fixation of the Ileo-sacral joint, and fusion of the foot and ankle.

For the 7.3mm Cannulated Cancellous Screw, indications for use are for fractures of the calcaneus, femoral neck, tibiaplateau, and of the sacrum and the articular cavity of the hip joint, fusion of the foot and ankle, fixation of the Ileo-sacral joint, and metaphyseal fractures of the distal femur and distal tibia.

The system(s) is not intended for spinal use.

BASES OF SUBSTANTIAL EQUIVALENCE:

The I.T.S. Straight Plate with Angular Stability is substantially equivalent to the Smith & Nephew, Zimmer, Synthes, Acumed, and I.T.S. GmbH stabilizing bone plate systems. The I.T.S. Screw System is substantially equivalent to the I.T.S. GmbH, Zimmer, and Synthes cannulated screw and washer systems.

SUMMARY OF SAFETY AND EFFECTIVENESS:

The I.T.S. Straight Plate with Angular Stability and Screw System is shown to be safe and effective for use in fracture fixation of small and long bones in the body.

6.5 mm Cannulated Cancellous Screw

Thread 22 mm



PROPERTIES

- back-tapping flank
- 6.5 mm
- length from 25 to 120 mm in 5 mm steps
- 3.2 mm guide wire

INDICATIONS

- fractures of the femoral neck
- metaphysal fractures of the distal femur and distal tibia
- fractures of the tibial plateau
- fractures of the sacrum and the articular cavity of the hip joint
- fixation of the Ileo-sacral joint
- fusion of foot and ankle

6.5 mm Cannulated SCFE Screw

Transcutaneous Screwing of the Slipped Capital Femoral Epiphysis



PROPERTIES

- 10 mm thread
- 3.2 mm guide wire
- Back-tapping flank
- Large countersink head for easy removableness (trocar 7.5 mm / 30°)



INDICATIONS

- loosening of the epiphysis femoralis

GAUGE 2.
WITH CL

7.3 mm Cannulated Cancellous Screw

fully threaded



PROPERTIES

- 7.3 mm
- length from 50 to 120 mm in 5 mm steps
- 3.2 mm guide wire

INDICATIONS

- fractures of the calcaneus
- fractures of the femoral neck

variable thread



PROPERTIES

- back-tapping flank
- 7.3 mm
- length from 40 to 120 mm in 5 mm steps
- 3.2 mm guide wire

INDICATIONS

- fractures of the femoral neck
- fractures of the tibial plateau
- fusion of foot and ankle
- fixation of Ileo-sacral joint
- fractures of the sacrum and the articular cavity of the hip joint
- metaphysal fractures of the distal femur and distal tibia

16 mm thread



PROPERTIES

- back-tapping flank
- 7.3 mm
- length from 80 to 120 mm in 5 mm steps
- 3.2 mm guide wire







DEPT




Accessories

FLAT WASHER

			
4,0		6,5	
36431 (flat)	36432 (curved)	36651 (flat)	36653 (flat)
∅ internal ∅ external	∅ internal ∅ external	∅ internal ∅ external	∅ internal ∅ external
4.5 11	4.5 11	7 14	7 12.6
for 4.0mm screws	for 4.0mm screws	for 6.5mm screws	for 6.5mm screws
	36652 (curved)	36731 (flat)	36732 (curved)
	∅ internal ∅ external	∅ internal ∅ external	∅ internal ∅ external
	7 14	7.5 15	7.5 14
	for 6.5mm screws	for 7.3mm screws	for 7.3mm screws

Washers are adjusted to the screw head and improve the formfitting of the screw head on the bone.

DEPTH GAUGE



Depth gauge for calibrated guide wire ø 3.2 mm

item no.

59322



Depth gauge for calibrated guide wire ø 1.6 mm

59162



Depth gauge for plate systems

59022



AAMI

Standards

and

Recommended

Practices

Volume 1.2:
Sterilization, Part 2—Hospital equipment and industrial process control

Documents developed under the auspices of the
AAMI Sterilization Standards Committee

by:

Biological Indicators Working Group
Dry Heat Sterilization Working Group
Ethylene Oxide Residuals Working Group
Hospital Ethylene Oxide Sterilizer Working Group
Hospital Steam Sterilizer Working Group
Industrial Ethylene Oxide Sterilization Working Group
Industrial Steam Sterilization Working Group
Microbiological Methods Working Group
Returned Devices Decontamination Working Group
Radiation Sterilization Working Group
Reusable Devices Resterilization Working Group

Association for the Advancement of Medical Instrumentation

Published by the
Association for the Advancement of Medical Instrumentation
3330 Washington Boulevard, Suite 400
Arlington, Virginia 22201

© 1995 by the
Association for the Advancement of Medical Instrumentation

All Rights Reserved

No part of this publication may be reproduced in any form, in an electronic retrieval system
or otherwise, without prior permission of the publisher.

Printed in the United States of America

ISBN 1-57020-035-1

170

Contents — Sterilization, Part 2

	DESIGNATION*	PAGE
Contents — Sterilization, Part 1		v
Foreword		vii
Introduction		viii
How to use this book		x
Safe handling of biologically contaminated medical devices in nonclinical and clinical settings (AAMI Technical Information Report)	TIR10	415
BIER/EO gas vessels, 2nd edition (American National Standard)	ST44	443
BIER/Steam vessels, 2nd edition (American National Standard)	ST45	461
Biological indicators for ethylene oxide sterilization processes in health care facilities (American National Standard)	ST21	475
Biological indicators for saturated steam sterilization processes in health care facilities (American National Standard)	ST19	491
Guideline for the use of ethylene oxide and steam biological indicators in industrial sterilization processes (American National Standard)	ST34	507
Hospital steam sterilizers, 3rd edition (American National Standard)	ST8	539
Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 2nd edition (American National Standard)	ST24	561
Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for device manufacturers (AAMI Technical Information Report)	TIR12	581
Medical devices—Validation and routine control of ethylene oxide sterilization, 3rd edition (American National Standard)	11135	643
Determining residual ethylene chlorohydrin and ethylene glycol in medical devices (American National Standard)	ST30	671

*This (alpha)numeric designation appears in the black tab on the side of each page of all documents.

Annex B

Sterilization cycles available in health care facilities

This annex describes sterilization cycles that are currently available for use in health care facilities or that are expected to be available for hospital applications by the end of 1994. Not all of the processes listed are in wide use at the time of publication of this report, and not all of them are suitable for any given medical device. It is important for medical device manufacturers to consult with appropriate sterilizer manufacturers for further information on the characteristics of new processes and their suitability for particular devices. It should be noted that pressure parameters are not given here for the various processes, some of which operate under significant positive or negative pressure. The processes and cycle parameters are listed for information only, and their inclusion here does not imply endorsement by AAMI.

Steam sterilization

Due to increasing cost constraints and, especially, to concerns about occupational exposure to toxic residuals of some gaseous and liquid chemical sterilization processes, sterilization by saturated steam is the preferred reprocessing method in health care facilities. Metal surgical instruments and other heat-stable medical devices are commonly sterilized by this method. Typical parameters for the generic cycles available for use in health care facilities are as follows:

Gravity-displacement steam sterilization

Wrapped items:

Temperature: 270° F to 275° F (132° C to 135° C)
Exposure time: 10 to 15 minutes

Temperature: 250° F to 254° F (121° C to 123° C)
Exposure time: 15 to 30 minutes

Unwrapped items ("flash" sterilization):

Temperature: 270° F (132° C)
Exposure time: 3 minutes (metal instruments only)
10 minutes (mixed porous, nonporous items)

Prevacuum steam sterilization

Wrapped items:

Temperature: 270° F to 275° F (132° C to 135° C)
Exposure time: 3 to 4 minutes

Unwrapped items ("flash" sterilization):

Temperature: 270° F (132° C)
Exposure time: 3 minutes (metal instruments only)
4 minutes (mixed porous, nonporous items)

TIR12

Steam-flush pressure-pulse steam sterilization

Wrapped items:

Temperature: 250° F to 254° F (121° C to 123° C)
Exposure time: 20 minutes

Unwrapped items:

Temperature: 270° F to 275° F (132° C to 135° C)
Exposure time: 3 to 4 minutes

Ethylene oxide sterilization

Ethylene oxide (EO) sterilization is commonly used to sterilize items that cannot withstand high temperatures, such as medical devices composed entirely or in part of plastic. Currently, health care facilities generally use either 100% ethylene oxide or a sterilant mixture consisting of 12% ethylene oxide and 88% chlorofluorocarbon-12 (CFC-12); CFC-12 is used as a diluent to reduce the flammability of ethylene oxide. Because of concerns about the harmful effects of CFCs on the ozone layer of the earth's atmosphere, alternative EO sterilant mixtures are under development and some have begun to be introduced in health care facilities. Such mixtures are not in wide use at present, but will become more prevalent in the future as regulatory constraints reduce the production and use of CFC-12. Typical sterilization cycle parameters for the commonly used 100% EO and 12/88 EO/CFC-12 methods are given below, along with the parameters for two of the newer mixtures:

100% EO

Concentration EO: 883 milligrams per liter (mg/L)
Temperature: 131° F (55° C) or 99° F (37° C)
Exposure time: 60 to 250 minutes (depending on temperature)
Humidity: 70% RH minimum

Concentration EO: 725 mg/L
Temperature: 131° F (55° C) or 99° F (37° C)
Exposure time: 60 to 180 minutes (depending on temperature)
Humidity: 70% RH minimum

510(k) Submission
I.T.S. Implantat-Technologie-Systeme GmbH
Pelvic Reconstruction System (PRS)

Section X: Labeling

See attached sample labeling for various components of the I.T.S. Pelvic Reconstruction System (PRS) designating the manufacturer and address, product catalog number, product description, quantity in package, material, lot/production number traceability, and sterilization designation as **NONSTERILE**.

A package insert for the I.T.S. Pelvic Reconstruction System (PRS) is currently in development that will disclose the Indications, Warnings, Precautions, Contraindications, Sterilization as NONSTERILE as labeled (with recommendations for in-hospital gravity autoclave per AAMI guidelines), and Materials. Also, the statement "Federal Law (USA) restricts this device to sale by or on order of a board certified physician" will be shown.

see
Section 11 RPS
R/22/2006

Sample Labeling for the I.T.S. PRS Sacral Rod




Implant-Technology-Systems
Aupal 28 A-8301 Laßnitzhöhe
Tel.net.: 0316/211-21-0 Fax 20
Tel. int.: +43/316/211-21-0 Fax 20

Made in Austria

Gewindestange threaded rod

SIZE L=200mm QTY 1
MAT TiAl6V4 LOT 65/6056 UNSTERIL NON-STERILE
REF 39401-01



CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a board certified physician.



Implant-Technology-Systems
Aupal 28 A-8301 Laßnitzhöhe
Tel.net.: 0316/211-21-0 Fax 20
Tel. int.: +43/316/211-21-0 Fax 20

Austria Mutter Nut

SIZE SW 10
MAT Ti 6Al 4V
REF 39401-03
QTY 1
LOT 65/6055



CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a board certified physician.



Implant-Technology-Systems
Aupal 28 A-8301 Laßnitzhöhe
Tel.net.: 0316/211-21-0 Fax 20
Tel. int.: +43/316/211-21-0 Fax 20

Austria

Beilagscheibe, keilförmig, gerippt Washer, wedge-shaped, ripped

MAT TiAl6V4
REF 39401-02
QTY 1
LOT 26/F14963



CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a board certified physician.



Implant-Technology-Systems
Aupal 28 A-8301 Laßnitzhöhe
Tel.net.: 0316/211-21-0 Fax 20
Tel. int.: +43/316/211-21-0 Fax 20

Austria

Sicherungsmutter Securing nut

SIZE SW 7
MAT Ti 6Al 4V ELI
REF 39401-04
QTY 1
LOT 65/6056



CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a board certified physician.

Sample Labeling for the I.T.S. PRS 7.3mm Cannulated SIJ Traction Screw & Washers



Implant-Technology-Systems
Aupal 28 A-8301 Laßnitzhöhe
Tel.nat.: 0316/211-21-0 Fax 20
Tel. int.: +43/316/211-21-0 Fax 20

Made in Austria

7,3mm kanülierte
Spongiosaschraube, variables Gewinde
7,3mm Cannulated
cancellous screw, variable thread

SIZE D=7,3mm; L=95mm

MAT Ti 6Al 4V ELI

REF 31734-95

QTY 1

LOT 26/F14877

UNSTERIL
NON-STERILE



CAUTION: Federal Law (USA) restricts this device to sale
by or on the order of a board certified physician.



Implant-Technology-Systems
Aupal 28 A-8301 Laßnitzhöhe
Tel.nat.: 0316/211-21-0 Fax 20
Tel. int.: +43/316/211-21-0 Fax 20

Made in Austria

Unterlegscheibe für
7,3mm Schraube, gewölbt
Cupped Washer for
7,3mm Screw

MAT Ti 6Al4V

REF 36732

QTY 1

LOT 65/6015

UNSTERIL NON-STERILE



CAUTION: Federal Law (USA) restricts this device to sale
by or on the order of a board certified physician.



Implant-Technology-Systems
Aupal 28 A-8301 Laßnitzhöhe
Tel.nat.: 0316/211-21-0 Fax 20
Tel. int.: +43/316/211-21-0 Fax 20

Made in Austria

Unterlegscheibe für
7,3mm Schraube, flach
Flat Washer for
7,3mm Screw

MAT Ti 6Al4V

REF 36731

QTY 1

LOT 65/6017

UNSTERIL NON-STERILE



CAUTION: Federal Law (USA) restricts this device to sale
by or on the order of a board certified physician.

Sample Labeling for the I.T.S. PRS Multiple Type - Pelvic Plating System



Implant-Technology-Systems
Austal 28 A-8301 Laßnitzhöhe
Tel.nat.: 0316/211-21-0 Fax 20
Tel. int.+43/316/211-21-0 Fax 20

Austria

Gerade Beckenrekonstruktionsplatte

Straight pelvise reconstruction plate

SIZE	12 Loch; 12 holes
MAT	Titan; CP Titanium Grade 2 per ASTM F67-00
REF	21181-12
QTY	1
LOT	26/F14965



CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a board certified physician.

License DE 4343117, EP 1143 867, US 6322562



Implant-Technology-Systems
Austal 28 A-8301 Laßnitzhöhe
Tel.nat.: 0316/211-21-0 Fax 20
Tel. int.+43/316/211-21-0 Fax 20

Austria

Gebogene Beckenrekonstruktionsplatte

Pelvise reconstruction plate

SIZE	8 Loch; 8 holes / R108
MAT	Titan; CP Titanium Grade 2 per ASTM F67-00
REF	21195-8
QTY	1
LOT	26/F14981



CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a board certified physician.

License DE 4343117, EP 1143 867, US 6322562



Implant-Technology-Systems
Austal 28 A-8301 Laßnitzhöhe
Tel.nat.: 0316/211-21-0 Fax 20
Tel. int.+43/316/211-21-0 Fax 20

Austria

Gebogene Beckenrekonstruktionsplatte

Pelvise reconstruction plate

SIZE	8 Loch; 8 holes / R88
MAT	Titan; CP Titanium Grade 2 per ASTM F67-00
REF	21194-8
QTY	1
LOT	26/F14972



CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a board certified physician.

License DE 4343117, EP 1143 867, US 6322562



Implant-Technology-Systems
Austal 28 A-8301 Laßnitzhöhe
Tel.nat.: 0316/211-21-0 Fax 20
Tel. int.+43/316/211-21-0 Fax 20

Austria

Symphysenplatte

Symphysis plate

SIZE	4 Loch; 4 holes
MAT	Titan; CP Titanium Grade 2 per ASTM F67-00
REF	21161-4
QTY	1
LOT	65/6079



CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a board certified physician.

License DE 4343117, EP 1143 867, US 6322562



Implant-Technology-Systems
Austal 28 A-8301 Laßnitzhöhe
Tel.nat.: 0316/211-21-0 Fax 20
Tel. int.+43/316/211-21-0 Fax 20

Austria

Symphysenplatte

Symphysis plate

SIZE	6 Loch; 6 holes
MAT	Titan; CP Titanium Grade 2 per ASTM F67-00
REF	21161-6
QTY	1
LOT	65/6079



CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a board certified physician.

License DE 4343117, EP 1143 867, US 6322562

Sample Labeling for the I.T.S. PRS Multiple Type - Pelvic Plating System



Implant-Technology-Systems
Austal 28 A-8301 Laßnitzhöhe
Tel.nat.: 0316/211-21-0 Fax 20
Tel. int.: +43/316/211-21-0 Fax 20

Austria

ISG Platte, geschlossen

SIJ-plate, closed

SIZE	4 Loch; 4 holes
MAT	Titan; CP Titanium Grade 2 per ASTM F67-00
REF	21173-4
QTY	1
LOT	26/F15056



CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a board certified physician.

License DE 4343117, EP 1143 867, US 6322562



Implant-Technology-Systems
Austal 28 A-8301 Laßnitzhöhe
Tel.nat.: 0316/211-21-0 Fax 20
Tel. int.: +43/316/211-21-0 Fax 20

Austria

ISG Platte, L-Form, links

SIJ-plate, L-shape, left

SIZE	5 Loch; 5 holes
MAT	Titan; CP Titanium Grade 2 per ASTM F67-00
REF	21172-5
QTY	1
LOT	26/F15054



CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a board certified physician.

License DE 4343117, EP 1143 867, US 6322562



Implant-Technology-Systems
Austal 28 A-8301 Laßnitzhöhe
Tel.nat.: 0316/211-21-0 Fax 20
Tel. int.: +43/316/211-21-0 Fax 20

Austria

4,5mm Kortikalisschraube

4,5mm Cortical Screw

SIZE	D=4,5mm; L=40mm
MAT	Ti 6Al 4V
REF	32455-40
QTY	1
LOT	26/F13834



CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a board certified physician.



Implant-Technology-Systems
Austal 28 A-8301 Laßnitzhöhe
Tel.nat.: 0316/211-21-0 Fax 20
Tel. int.: +43/316/211-21-0 Fax 20

Austria

**Spongiaschraube für
Beckenrekonstruktionsplatten**

**Anglestable cancellous screw
for pelvise reconstruction plate**

SIZE	D=5,9mm; L=40mm
MAT	Ti 6Al 4V ELI
REF	30591-40
QTY	1
LOT	26/F15309



CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a board certified physician.

License DE 4343117, EP 1143 867, US 6322562



Implant-Technology-Systems
Austal 28 A-8301 Laßnitzhöhe
Tel.nat.: 0316/211-21-0 Fax 20
Tel. int.: +43/316/211-21-0 Fax 20

Austria

**Spongiosa-
Stabilisierungsschraube**

**Cancellous-
stabilization screw**

SIZE	D=5,9mm; L=40mm
MAT	Ti 6Al 4V ELI
REF	37592-40
QTY	1
LOT	26/F15322



CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a board certified physician.

License DE 4343117, EP 1143 867, US 6322562



Indications for Use

510(k) NUMBER: _____

DEVICE NAME: **PELVIC RECONSTRUCTION SYSTEM (PRS)**

The ***intended use*** of the I.T.S. Pelvic Reconstruction System (PRS) is to stabilize one or more pelvic bone fractures in the 'pelvic ring' area of the pelvis of an adult or pediatric patient which include the use of the following pelvic fracture fixation systems.

Indications for use of the I.T.S. PRS Low Profile - Multiple Type – Pelvic Plating System include:

- 1) Fracture reconstruction of the acetabulum, sacrum, ilium, and entire pelvic ring,
- 2) Revision surgery of pseudarthroses, non-unions and mal-unions,
- 3) Ilio-Iliac distance osteosynthesis,
- 4) Osteotomies,
- 5) Arthrodesis,
- 6) Sacroiliac joint dislocations, and
- 7) Symphysis pubis disruptions

Indications for use of the I.T.S. PRS Sacral Rod System include:

- 1) Fixation of fractures of the posterior pelvis,
- 2) Fixation of fractures of the posterior iliac spine,
- 3) Fixation of fractures of the posterior inferior iliac spine,
- 4) Dorsal stabilization of the posterior pelvic ring for unstable pelvic ring injuries,
- 5) Fixation of sacral fractures, and
- 6) Fracture dislocations of the sacro-iliac joint

Indications for use of the I.T.S. PRS 7.3mm Cannulated SIJ Traction Screw & Washer include fracture fixation of pelvic bone fractures where indicated and for periacetabular osteotomies.

has not been studied in

The system(s) ~~is not intended for spinal use,~~ and is not intended for use in vertebral column fracture or fusion procedures.

RST 12/22/2006

Prescription Use AND/OR Over-The-Counter-Use _____

(Per 21 CFR 801 Subpart D)

(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

179

Section XII: 510(k) Summary of Safety and Effectiveness

SAFE MEDICAL DEVICES ACT OF 1990

510(k) Summary

NAME OF FIRM: I.T.S. Implantat-Technologie-Systeme GmbH.
Autal 28.
Lassnitzhoehe A – 8301
AUSTRIA

510(k) FIRM CONTACT: Al Lippincott
Engineering Consulting Services, Inc.
3150 E. 200th St.
Prior Lake, MN 55372

TRADE NAME: Pelvic Reconstruction System (PRS)

COMMON NAME: Pelvic Bone Fixation Set

CLASSIFICATION: Plate, Fixation, Bone (see 21 CFR, Sec. 888.3030),
Pin, Fixation, Threaded (see 21 CFR, Sec. 888.3040)
Screw, Fixation, Bone (see 21 CFR, Sec. 888.3040),
Washer, Bolt, Nut (see 21 CFR, Sec. 888.3030)

DEVICE PRODUCT CODE: HRS

SUBSEQUENT PRODUCT CODE: JDW, HWC, HTN

SUBSTANTIALLY EQUIVALENT DEVICES Synthes 3.5mm Low Profile Pelvic Reconstruction (**K031573**)
Stryker Trauma Pelvic Set (**K001614**)
Zimmer Pelvic Reconstruction Set
Synthes Sacral Bar System (**K001720**)
I.T.S. GmbH FR.O.H. Calcaneus Repair System (**K051642**)
Synthes 3.9mm Pelvic Screw (**K013044**)
I.T.S. Straight Plate with Angular Stability & Screw System (**K060156**)

DEVICE DESCRIPTION: The I.T.S. Pelvic Reconstruction System (PRS) encompasses a number of fracture fixation subsystems (multiple pelvic plate designs, sacral threaded rod, and cannulated screw & washer) for fracture fixation and reconstruction of pelvic ring fractures in the pelvis.

The I.T.S. PRS Low Profile – Multiple Type – Pelvic Plating System consists of the following plate types: 1) A Straight Plate at a 4.0mm thickness with 10 to 14 hole length sizes, 2) A Curved Plate with a 108mm radius at a 2.5mm thickness with 4 to 16 hole length sizes, 3) A Curved Plate with a 88mm radius at a 2.5mm thickness also with 4 to 16 hole length sizes, 4) A Symphysis Plate at a 4.0mm thickness in both a 4 & 6 hole size, 5) A Sacroiliac-Joint (SIJ) L-Shaped Plate at a 2.5mm thickness in a left and right 5 hole size, and 6) A SIJ Closed Plate at a 2.5mm thickness in a 4 hole size.

All plate designs are low profile in thickness and made from CP Titanium. The PRS Pelvic Plating System also encompasses a number

of Cancellous (5.9mm std.compression & 5.9mm Locking screws) and Cortical (4.5mm std.compression screw) screw types and length sizes. All bone screws are pre-drilling and self-tapping in design and manufactured from high strength 6-4 Alloyed Titanium. All components (plates & screws) have a TIODIZE II surface treatment preparation.

The I.T.S. PRS Sacral Rod System consists of a threaded pin, wedge-shaped washers, and nut/locknut design with one end of the threaded pin having a trocar point. The assembled unit is used to stabilize ilio-iliac posterior pelvic ring disruption injuries. All components are manufactured from high strength 6-4 Alloyed Titanium material and have the TIODIZE II surface preparation.

The I.T.S. PRS 7.3mm Cannulated SIJ Traction Screw & Washer is a pre-drilling, self-tapping, and back-tapping screw design for guided reduction of pelvic bone fractures. Washers (flat & curved) are available for use with the screw. All 7.3mm Cannulated Screws & washers are manufactured from high strength 6-4 Alloyed Titanium and have the TIODIZE II surface preparation.

INTENDED USE:

The *intended use* of the I.T.S. Pelvic Reconstruction System (PRS) is to stabilize one or more pelvic bone fractures in the pelvic ring area of the pelvis of an adult or pediatric patient which include the use of the following pelvic fracture fixation systems.

Indications for use of the I.T.S. PRS Low Profile - Multiple Type – Pelvic Plating System include:

- 1) Fracture reconstruction of the acetabulum, sacrum, ilium, and entire pelvic ring,
- 2) Revision surgery of pseudarthroses, non-unions and mal-unions,
- 3) Ilio-Iliac distance osteosynthesis,
- 4) Osteotomies,
- 5) Arthrodesis,
- 6) Sacroiliac joint dislocations, and
- 7) Symphysis pubis disruptions

Indications for use of the I.T.S. PRS Sacral Rod System include:

- 1) Fixation of fractures of the posterior pelvis,
- 2) Fixation of fractures of the posterior iliac spine,
- 3) Fixation of fractures of the posterior inferior iliac spine,
- 4) Dorsal stabilization of the posterior pelvic ring for unstable pelvic ring injuries,
- 5) Fixation of sacral fractures, and
- 6) Fracture dislocations of the sacro-iliac joint

Indications for use of the I.T.S. PRS 7.3mm Cannulated SIJ Traction Screw & Washer include fracture fixation of pelvic bone fractures where indicated and for periacetabular osteotomies.

The system(s) is not intended for spinal use. [△] see Section 11.

RAT 02/24/2006

BASIS OF SUBSTANTIAL EQUIVALENCE:

The I.T.S. Pelvic Reconstruction System (PRS) is substantially equivalent to the Synthes, Zimmer, and Stryker pelvic reconstruction systems and the I.T.S. GmbH stabilizing bone plating systems.

SUMMARY OF SAFETY AND EFFECTIVENESS:

The I.T.S. Pelvic Reconstruction System (PRS) is shown to be safe and effective for use in 'pelvic ring' bone fracture fixation of the pelvis.

JUN - 6 2003

K031573

3.0 510(k) Summary

Sponsor: Synthes (USA)
1690 Russell Road
Paoli, PA 19301
(610) 647-9700

Contact: Bonnie Smith

Device Name: Synthes 3.5 mm Low Profile Pelvic Reconstruction Plate

Classification: Class II, § 888.3030 – Single/multiple component metallic bone fixation appliances and accessories

Predicate Device: Synthes 3.5 mm LCP Reconstruction Plate

Device Description: Synthes 3.5 mm Low Profile Pelvic Reconstruction Plates consist of curved plates and straight plates that are contourable in three planes. The plates have a smooth, rounded profile. Plate holes are slightly oval, allowing 30 degrees of angulation in all directions. The 3.5 mm Low Profile Reconstruction Plate accepts Synthes 3.5 mm Cortex Screws and Synthes 3.5 mm Pelvic Screws.

Intended Use: Synthes 3.5 Low Profile Pelvic Reconstruction Plate is intended for pelvic and acetabular reconstructive surgery and fracture fixation of the distal humerus, clavicle, and scapula.

Materials: Stainless steel


00005

3.5 mm Low Profile Pelvic System

The 3.5 mm Low Profile Pelvic System provides additional features and benefits to the pelvic and acetabular reconstructive surgeon.

Plate Features

- Low-profile plate improves contourability for easy intraoperative bending, when necessary
- Improved screw angulation
- Low-profile plate and edges to minimize soft tissue irritation
- Variety of precontoured shapes to minimize intraoperative contouring
- Compatible with both existing and new instruments and screws
- Implant quality 316L stainless steel
- Coaxial Combi holes in locking reconstruction straight and J-plates, and in Pubic Symphysis Plates, allow use of locking screws or cortex screws in the same round conical plate hole



3.5 mm Low Profile Pelvic System Instruments

- New instruments specifically designed to provide improved function with the 3.5 mm Low Profile Reconstruction Plates.

Indications

The Synthes 3.5 mm Low Profile Reconstruction Plate is intended for pelvic and acetabular reconstructive surgery and fracture fixation of the distal humerus, clavicle, and scapula.

- Designed for easy contouring
- Edge designed to minimize soft tissue irritation
- 316L stainless steel

3.5 mm Low Profile Reconstruction Plates (straight)

- 3–20 holes
- Improved screw angulation
- Available with or without coaxial Combi holes

3.5 mm Low Profile Curved Reconstruction Plates

- 88 mm or 108 mm radius of curvature
- 6–16 holes
- Improved screw angulation

3.5 mm Wide Angle Low Profile Reconstruction Plates (straight)

- Larger diameter hole in thicker plate allows maximum screw angulation
- 3–20 holes

3.5 mm Low Profile Reconstruction J-Plates

- Precontoured to fit the pelvic brim
- 88 mm radius of curvature in curved portion of plate
- Left and right plates
- 10–16 holes
- Available with or without coaxial Combi holes


3.5 mm Pubic Symphysis Plates

- 4 and 6 coaxial Combi holes
- Precontoured to fit the pubic symphysis
- 60 mm radius of curvature
- 3 suture holes
- Cold-worked for improved strength

Screws

3.5 mm Cortex Screws and 3.5 mm Pelvic Cortex Screws

- Deep hexagonal recess minimizes stripping of the screwhead or driver
- 10 mm to 150 mm lengths



3.5 mm Low Profile Pelvic System Implant and instrument configurations may be customized. The customer may choose which of the following trays to combine.

Trays can be combined with any one of the graphic cases or frame case. This is listed below:

690.411	Single Level Graphic Case Base, 510 mm x 250 mm
690.412	Two Level Graphic Case Base, 510 mm x 250 mm
690.413	Three Level Graphic Case Base, 510 mm x 250 mm
690.422	Lid for 3.5 mm Low Profile Pelvic Retractors
690.426	Lid for 3.5 mm Low Profile Pelvic System Implant Set
690.427	Lid for 3.5 mm Low Profile Pelvic System Instruments
690.428	Lid for 3.5 mm Low Profile Pelvic System Reduction Instruments
690.429	Lid for 3.5 mm Low Profile Pelvic System
690.421	Lower Tray for 3.5 mm Low Profile Pelvic Retractors
690.423	Tray for 3.5 mm Low Profile Reconstruction J-Plates
690.424	Tray for 3.5 mm Low Profile Wide Angle Reconstruction Plates
690.431	Tray for 3.5 mm Low Profile Pelvic System Instruments
690.432	Upper Tray for 3.5 mm Low Profile Pelvic Retractors
690.433	Lower Tray for 3.5 mm Low Profile Pelvic System Reduction Instruments
690.434	Upper Tray for 3.5 mm Low Profile Pelvic System Reduction Instruments
690.435	Tray for 3.5 mm Locking Low Profile Reconstruction Plates
690.436	Tray for 3.5 mm Low Profile Reconstruction Plates
690.437	Tray for 3.5 mm Low Profile Pelvic System Cortex Screws
690.438	Screw Rack for 3.5 mm Low Profile Locking Screws
690.439	Tray for 3.5 mm Low Profile Pelvic System StarDrive Cortex Screws
690.916	Full Tray for 510 mm x 250 mm graphic case
690.917	Full Tray Finger Mat for 510 mm x 250 mm graphic case

3.5 mm Low Profile Pelvic System Instrument Configurations (including cases and lids)

01.100.009	3.5 mm Low Profile Pelvic System Implant Set
01.100.008	3.5 mm Low Profile Pelvic System Reduction Instrument Set
01.100.004	3.5 mm Low Profile Pelvic System Retractor Set
01.100.013	3.5 mm Low Profile Pelvic System Instrument Set

3.5 mm Low Profile Pelvic System Tray Configurations

01.100.012	3.5 mm Low Profile Pelvic Reconstruction Plate Tray
01.100.027	3.5 mm Wide Angle Low Profile Reconstruction Plate Tray
01.100.031	3.5 mm Low Profile Reconstruction J-Plate Tray

3.5 mm Low Profile Pelvic System Tray Configurations

01.100.042	3.5 mm Locking Low Profile Reconstruction Plate Tray
01.100.112	3.5 mm Low Profile Pelvic System Screw Tray with Locking Screws
01.100.122	3.5 mm Low Profile Pelvic System Screw Rack for Locking Screws
01.100.132	3.5 mm Low Profile Pelvic System StarDrive Screw Tray without Locking Screws
01.100.142	3.5 mm Low Profile Pelvic System Screw Tray without Locking Screws

Instruments

03.100.011	Spoon Retractor, medium, 268 mm
03.100.012	Spoon Retractor, large, 323 mm
03.100.013	Sciatic Nerve Retractor
03.100.014	Sciatic Nerve Retractor, long
03.100.015	Malleable Retractor, 20 mm width
03.100.016	Malleable Retractor, 30 mm width
03.100.017	Malleable Retractor, 40 mm width
03.100.018	Straight Ball Spike, 6.5 mm ball diameter, 320 mm
03.100.019	Straight Ball Spike, long, 6.5 mm ball diameter, 420 mm
03.100.020	Straight Reduction Forceps, with pointed ball tips, 6.5 mm ball diameter, 257 mm
03.100.021	Angled Reduction Forceps, small, with pointed ball tips, 6.5 mm ball diameter, 195 mm
03.100.022	Angled Reduction Forceps, large, with pointed ball tips, 6.5 mm ball diameter, 240 mm
03.100.023	Reduction Forceps, with pointed ball tips, 6.5 mm ball diameter, 400 mm
03.100.024	Asymmetric Reduction Forceps, with pointed ball tips, 6.5 mm ball diameter, 404 mm
03.100.025	Low Profile Pelvic System Reduction Forceps, used with 3.5 mm screws, 250 mm
03.100.027	Round Disk, 6.5 mm hole diameter
03.100.028	Rectangular Disk, 6.5 mm hole diameter
03.100.029	Medium Bone Hook, long, 330 mm
03.100.030	Large Bone Hook, long, 330 mm
03.100.031	Bending Pliers, for 3.5 mm Low Profile Reconstruction Plates
03.100.032	Ratcheting Handle with quick coupling
03.100.033	3.5 mm Screwdriver Shaft, hexagonal, long
03.100.040	2.8 mm Threaded Drill Guide, short
03.100.045	3.5 mm StarDrive Screwdriver Shaft, T15, 250 mm long
	Bending Templates
03.100.034	for 3.5 mm Low Profile Reconstruction Plates, 20 holes
03.100.035	for 3.5 mm Low Profile Reconstruction J-Plates (88 mm radius), 16 holes
03.100.036	for 3.5 mm Low Profile Curved Reconstruction Plates, 88 mm radius, 16 holes
03.100.037	for 3.5 mm Low Profile Curved Reconstruction Plates, 108 mm radius, 16 holes

3.5 mm Low Profile Pelvic Plates

Nontlocking		Locking		Locking	
Part No.	Part No.	Part No.	Part No.	Part No.	Part No.
245.023	02.100.07	12 holes	120 mm	02.100.07	12 holes
245.024	02.100.08	14 holes	120 mm	02.100.08	14 holes
245.025	02.100.09	16 holes	120 mm	02.100.09	16 holes
245.026	02.100.09	18 holes	120 mm	02.100.09	18 holes
245.027	02.100.10	18 holes	130 mm	02.100.10	18 holes
245.028	02.100.10	20 holes	130 mm	02.100.10	20 holes
245.029	02.100.10	22 holes	130 mm	02.100.10	22 holes
245.030	02.100.11	22 holes	140 mm	02.100.11	22 holes
245.031	02.100.11	24 holes	140 mm	02.100.11	24 holes
245.032	02.100.11	26 holes	140 mm	02.100.11	26 holes
245.033	02.100.11	28 holes	140 mm	02.100.11	28 holes
245.034	02.100.11	30 holes	140 mm	02.100.11	30 holes
245.035	02.100.11	32 holes	140 mm	02.100.11	32 holes
245.036	02.100.11	34 holes	140 mm	02.100.11	34 holes
245.038	02.100.11	36 holes	140 mm	02.100.11	36 holes
245.039	02.100.12	20 holes	160 mm	02.100.12	20 holes

3.5 mm Low Profile Curved Reconstruction Plates

88 mm Radius	108 mm Radius	Holes	Length	Part No.	Holes	Part No.
245.906	245.876	6 holes	40 mm	245.906	6 holes	
245.908	245.878	8 holes	40 mm	245.908	8 holes	
245.910	245.880	10 holes	40 mm	245.910	10 holes	
245.912	245.882	12 holes	40 mm	245.912	12 holes	
245.914	245.884	14 holes	40 mm	245.914	14 holes	
245.916	245.886	16 holes	40 mm	245.916	16 holes	

Screws

Part No.	Description	Length	Part No.	Holes	Part No.
204.640-204.750	3.5 mm Pelvic Cortical self-tapping, 40°	10-50 mm	02.100.36	12 holes	02.100.36
204.810-204.838	3.5 mm Cortex self-tapping, 10°	10-30 mm	02.100.36	14 holes	02.100.36
			02.100.36	16 holes	02.100.36

219.98 Washer, 7.0 mm

*5 mm increments
**2 mm increments

Part No.	Description	Length	Part No.	Holes	Part No.
02.100.00	3.5 mm Pelvic Cortical self-tapping, 40°	10-50 mm	02.100.00	12 holes	02.100.00
02.100.00	3.5 mm Cortex self-tapping, 10°	10-30 mm	02.100.00	14 holes	02.100.00

SYNTHES (USA)
1302 Wrights Lane East
West Chester, PA 19380
Telephone: (610) 719-5000
To order: (800) 523-0322
Fax: (610) 251-9056

SYNTHES (CANADA)
2800 Meadowcroft
Mississauga, Ont.
Telephone: (905) 276-8000
Toll free: (800) 387-2222
Fax: (905) 567-4100





Inventory Control Form

3.5 mm Low Profile Pelvic System

Patient Data - For Charge Use

Name: _____

Date: _____

SYNTHES (USA)
To Order: (800) 523-0322

SYNTHES (Canada) Ltd.
To Order: (800) 668-1119

Screws

3.5 mm Pelvic Cortex Screws, self-tapping

LENGTH
204.640
204.645
204.650
204.655
204.660
204.665
204.670
204.675
204.680
204.685
204.690
204.695
204.700
204.705
204.710
204.715
204.720
204.725
204.730
204.735
204.740
204.745
204.750

3.5 mm Locking Screws, self-tapping, with StarDrive recess

LENGTH
212.101
212.102
212.103
212.104
212.105
212.106
212.107
212.108
212.109
212.110
212.111
212.112
212.113
212.115
212.116
212.117
212.119
212.121
212.123
212.124
212.125
212.126
212.127
212.128
212.129
212.130
212.131

3.5 mm Cortex Screws, self-tapping, with T-15 StarDrive recess

LENGTH
02.200.010
02.200.012
02.200.014
02.200.016
02.200.018
02.200.020
02.200.022
02.200.024
02.200.026
02.200.028
02.200.030
02.200.032
02.200.034
02.200.036
02.200.038
02.200.040
02.200.042
02.200.044
02.200.045
02.200.046
02.200.048
02.200.050
02.200.055
02.200.060
02.200.065
02.200.070
02.200.075
02.200.080
02.200.085
02.200.090
02.200.095
02.200.100
02.200.105
02.200.110
02.200.115
02.200.120
02.200.125
02.200.130
02.200.135
02.200.140
02.200.145
02.200.150

3.5 mm Cortex Screws, self-tapping

LENGTH
204.810
204.812
204.814
204.816
204.818
204.820
204.822
204.824
204.826
204.828
204.830
204.832
204.834
204.836
204.838

3.5 mm Threaded Plug, StarDrive

(for protection of the 3.5 mm Threaded Plate holes during contouring)

02.200.003

188

Plates

3.5 mm Low Profile Reconstruction Plates

HOLES	LENGTH
3	39 mm
4	52 mm
5	65 mm
6	78 mm
7	91 mm
8	104 mm
9	117 mm
10	130 mm
11	143 mm
12	156 mm
13	169 mm
14	182 mm
15	195 mm
16	208 mm
18	234 mm
20	260 mm

3.5 mm Low Profile Reconstruction J-Plates (88 mm radius)

LEFT	RIGHT	HOLES	LENGTH
245.930	245.920	10	130 mm
245.932	245.922	12	156 mm
245.934	245.924	14	182 mm
245.936	245.926	16	208 mm

3.5 mm Low Profile Curved Reconstruction Plates

88 mm radius	108 mm radius	HOLES	LENGTH
245.906	245.876	6	78 mm
245.908	245.878	8	104 mm
245.910	245.880	10	130 mm
245.912	245.882	12	156 mm
245.914	245.884	14	182 mm
245.916	245.886	16	208 mm

3.5 mm Wide Angle Low Profile Reconstruction Plates

HOLES	LENGTH
3	39 mm
4	52 mm
5	65 mm
6	78 mm
7	91 mm
8	104 mm
9	117 mm
10	130 mm
11	143 mm
12	156 mm
13	169 mm
14	182 mm
15	195 mm
16	208 mm
18	234 mm
20	260 mm

Miscellaneous Implants

219.98	Washer, 7 mm
292.20	2.0 mm Kirschner Wire, 150 mm, 10 pack
294.68	6.0 mm Schanz Screw, spade point, 190 mm

Miscellaneous Instruments

Drill Bits, quick coupling

DIAMETER	LENGTH
310.288	2.8 mm
310.37	3.5 mm
315.92	2.5 mm
324.210	2.5 mm
324.211	2.8 mm

Tap

311.31	for 3.5 mm Cortex Screws, 175 mm, 140 mm calibration
--------	--

Bending Templates

03.100.034	for 3.5 mm, Low Profile Reconstruction Plates, 20 holes, 260 mm
03.100.035	for 3.5 mm, Low Profile Reconstruction J-Pates, 16 holes for 3.5 mm, Low Profile Curved Reconstruction Plates, 88 mm radius, 16 holes
03.100.036	for 3.5 mm, Low Profile Curved Reconstruction Plates, 108 mm radius, 16 holes
03.100.037	for 3.5 mm, Low Profile Curved Reconstruction Plates, 108 mm radius, 16 holes

3.5 mm Locking Low Profile Reconstruction J-Plates (88 mm radius)

LEFT	RIGHT	HOLES	LENGTH
02.100.361	02.100.360	10	130 mm
02.100.363	02.100.362	12	156 mm
02.100.365	02.100.364	14	182 mm
02.100.367	02.100.366	16	208 mm

3.5 mm Locking Low Profile Pubic Symphysis Plate with Coaxial Combi Holes

HOLES	LENGTH
4	57 mm
6	78 mm

AUG 4 2000

510(k) PREMARKET NOTIFICATION
SUMMARY OF SAFETY AND EFFECTIVENESS
Stryker Trauma Pelvic Set

K001614

Submission Information

**Name and Address of the Sponsor
of the 510(k) Submission:**

Howmedica Osteonics Corp.
59 Route 17
Allendale, NJ 07401-1677
201-825-4900

Contact Person:

Mary-Catherine Dillon
Regulatory Affairs Team Member

Date Summary Prepared:

March 23, 2000

Device Identification

Proprietary Name:

Stryker Trauma Pelvic Set

Common Name:

Pelvic Set

Classification Name and Reference:

Plate, Fixation, Bone
21 CFR §888.3030

Predicate Device Identification

The Stryker Trauma Pelvic Set is substantially equivalent to the Synthes Pelvic Implant Set.

Device Description

The Stryker Trauma Pelvic Set consists of 88mm radius and 108mm radius curved plates, straight pelvic plates, straight acetabulum plates, and symphysis-pubis plates. All curves and straight plate components are available in 10.5mm widths and 2.5mm thicknesses. The symphysis-pubis plates are also available in a 12.5mm width and 3.2mm thickness. The subject components vary in length from 22.5mm to 474.5mm. The system also includes 3.5mm diameter and 4.5mm diameter screws. All devices in the system are provided both sterile and non-sterile.

Intended Use

The Stryker Trauma Pelvic Set is indicated for:

- Fractures of the acetabulum, sacrum, illium, and entire pelvic ring
- Revision surgery of pseoduarthroses, non-unions and mal-unions

190


- Osteonmics
- Arthrodeses
- Sacroilic joint dislocations
- Symphysis pubis disruptures

Materials

The subject components are manufactured from stainless steel which conforms to ASTM F-138.

Statement of Technological Comparison

The subject components of the Stryker Trauma Pelvic Set are substantially equivalent in design and intended use to the predicate device offered by Synthes.



Trauma

External Fixation
 Hoffmann® II
 Hoffmann® II Micro
 Hoffmann® II Hybrid
 Hoffmann® II Compact™
 DJD® II
 Apex Pin Fixation
 Manotube® Triax™

Immediate Care
 Products
 Surgical Techniques

Hip Fracture/Int. Fixation
 Omega Plus™
 Anis III™
 Matta Pelvic System

IM Nails
 Gamma™ Locking Nail
 T2™ Tibia
 T2™ Femur
 T2™ Humerus

About Trauma
 Driving Directions
 Mailing Address

Find A Sales Rep

City of Fracture Internal Fixation

Stryker Trauma > Hip Fracture / Internal Fixation

Matta Pelvic Plating System

Pelvic surgery represents one of the most technically demanding specialties within orthopaedic trauma. The Matta Pelvic Plating System offers the highly specialized surgeon a complete and specific system for addressing fractures of the pelvis and acetabulum.

Joel Matta, MD, world renowned as a premiere pelvic surgeon, is the designer of this system. The system features specific refinements including some of the familiar concepts he learned as a Letournel Fellow.

The system also offers a comprehensive, single set of instrumentation including reductions forceps, clamps and sciatic nerve retractors.

INDICATIONS | SYSTEM CONFIGURATION | SYSTEM



INDICATIONS

- Fractures of the:
 - Acetabulum
 - Sacrum
 - Ilium
 - Pelvic Ring
- Revision Surgery of pseudoarthroses, non-unions and mal-unions
- Osteotomies
- Arthrodeses
- Sacroiliac Joint Dislocations
- Symphysis Pubis Disruptures

[Back to Top]

SYSTEM CONFIGURATION

- Material
 - Stainless Steel (Cold-Forged and Annealed)
 - Conforms to ISO/ASTM Standards: ISO 5832-1/ASTM F138
- Plating Options
 - Straight and curved plating options
 - Pre-curved Pelvic Plates (R108 and R88) to accommodate male and female anatomy
 - Low Screw Head Profile (Stainless Steel)
 - Fully Threaded, Self Tapping Cortical Screws;
Diameter: 3.5mm and 4.5mm
 - Cancellous Screw;
Diameter: 6.5mm; Lengths: 16mm and 32mm

[Back to Top]

192

SYSTEM HIGHLIGHTS

- Dedicated Symphysis Pubis Plate
- Annealed Acetabular Plate offers easier three dimensional bending
- Speciality Reduction Forceps and optimal Bone Clamp design offers improved visualization of the wound site
- Dedicated set of Instrumentation consists of Basic Instrument Tray, Reduction Instrument Tray, Implant Tray and Screw Rack

[Back to Top]



Symphysis Pub

193




Zimmer® Plate and Screw System

A New Focus on Trauma Solutions

The Zimmer Plate and Screw System offers complete implant and instrument portfolios for mini, small basic fragment and pelvic reconstruction systems, in sterilization cases featuring unequalled convenience, storage and flexibility.


MINI FRAGMENT SET

- Plates and screws in 1.5mm, 2.0mm, and 2.7mm sizes
- Extensive screw line with hex and cruciform head options, and standard or self-tapping flute options
- Sterilization case stores optional implants, eliminating the need to sterilize individual products




SMALL FRAGMENT SET

- Plates and screws in 3.5mm and 4.0mm sizes
- Full line of Compression Plates, 1/3 tubular plates, T plates, and Cloverleaf plates available
- Sterilization system combines flexibility and convenience with modular design and optional storage



BASIC FRAGMENT SET


- Plates and screws in 4.5mm and 6.5mm sizes; additional plate sizes, including 4.5mm reconstruction plates, can be stored in sterilization case base
- Sterilization container offers secure stacking to conserve table space



PELVIC RECONSTRUCTION SET

- Plates and screws in 3.5mm, 4.5mm, and 6.5mm sizes

- Comprehensive scope of implants included in sterilization case improves hospital efficiency by reducing the need to open additional sets
- Provides complete line of Reduction Clamps for pelvic fracture manipulations



[View Brochure](#)

Please refer to the package inserts for complete product information, including contraindications, warnings, precautions, and adverse effects.

Warning: This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.

Last updated 4/20/06

www.zimmer.com

© 2005 Zimmer, Inc.

Although accessible by others, the product information on the Zimmer site is intended for use by U.S. residents only unless designated by a national flag or regional symbol. This site may contain information related to various health, medical, and fitness conditions and their treatment. Such information is not meant to be a substitute for the advice provided by a physician or other medical professional. You should always consult a professional medical adviser, such as a physician, to determine what courses of treatment, if any, may be appropriate for you.

195

ZIMMER PLATES, SCREWS, AND INSTRUMENTS

- Plates and screws in 3.5mm and 4.0mm sizes
- A full line of Compression Plates, 1.3 Tubular plates, T-plates, and Cloverleaf plates is available
- New sterilization system combines flexibility and convenience with modular design and optional storage




SMALL FRAGMENT SET



PELVIC RECONSTRUCTION SET

- Plates and screws in 3.5mm, 4.5mm, and 6.5mm sizes. Comprehensive scope of implants included in the sterilization case improves hospital efficiency by reducing the need to open additional sets.
- Provides a complete line of Reduction Clamps for manipulation of pelvic fractures.



↑ Zimmer Pelvic Set



PELVIC RECONSTRUCTION SET

- Plates and screws in 3.5mm, 4.5mm, and 6.5mm sizes. Comprehensive scope of implants included in the sterilization case improves hospital efficiency by reducing the need to open additional sets.
- Provides a complete line of Reduction Clamps for manipulation of pelvic fractures.

AUG 24 2000

3.0 Summary of Safety and Effectiveness Information

SPONSOR: Synthes (USA)
1690 Russell Road
Paoli, PA 19301
(610) 647-9700
Contact: Angela Silvestri

DEVICE NAME: Synthes Sacral Bar System

CLASSIFICATION: 21 CFR 888.3040
Smooth or threaded metallic bone fixation fastener.

PREDICATE DEVICES: Zimmer Threaded Sacral Rod
Synthes Threaded Bolt


DEVICE DESCRIPTION: The Synthes Sacral Bar System consists of a threaded bar, washers, and nuts. The bars are fully threaded. One end of the bar has a trocar point to guide the bar through pre-drilled holes. The bars are available in lengths ranging from 120 to 260 mm, in 10 mm increments. The washers that are used with this system are oval shaped and are designed to slide freely along the bars. Both rounded and straight nuts are provided with this system; the rounded nuts mate with the washers to create compression, while the straight nuts are then added to wedge against the rounded nuts to maintain compression.

INTENDED USE: The Synthes Sacral Bar System is intended for fixation of fractures of the posterior pelvis, in areas of the posterior superior iliac spine and posterior inferior iliac spine, for sacral fractures and fracture-dislocations of the sacro-iliac joint.

MATERIAL: Implant quality stainless steel.

Sacral Bar Kit

TECHNIQUE GUIDE



Product Information

6.0 mm Threaded Sacral Bar Kits[◇]

199.931	160 mm
199.932	180 mm
199.933	200 mm
199.934	220 mm
199.935	240 mm
199.936	260 mm

Note:

Construct described in this Technique Guide requires two Sacral Bar Kits.

Each kit contains:

296.76x	6.0 mm Threaded Sacral Bar, 1 ea.
296.757	6.0 mm Sacral Nuts, rounded, 2 ea.*
296.758	6.0 mm Sacral Washers, 2 ea.*
296.759	6.0 mm Hexagonal Nuts, 2 ea.*
310.60J	6.0 mm Drill Bit, Jacobs chuck, 195 mm*

6.0 mm Threaded Sacral Bars*

296.761	160 mm
296.762	180 mm
296.763	200 mm
296.764	220 mm
296.765	240 mm
296.766	260 mm

Additionally Available

150.16	Compact Air Drive II Set or
105.957	Power Drive Set
321.16	Combination Wrench, 11 mm width across flats
321.20	Ratchet Wrench, 11 mm width across flats
388.72	Rod Cutter (<i>necessary to complete procedure</i>)
393.10	Universal Chuck with T-Handle
398.86	Pelvic Reduction Forceps, with pointed-ball tips, 400 mm length, extra large, long jaws, speed lock
511.20	Oscillating Drill Attachment

[◇] Kits available nonsterile or sterile-packed. Add "S" to catalog number to order sterile Sacral Bar Kit.

* Individual components may be purchased separately, nonsterile only.

K051642

AUG 22 2005

510(k) Summary of Safety and Effectiveness

SAFE MEDICAL DEVICES ACT OF 1990

510(k) Summary

NAME OF FIRM: I.T.S. Implantat-Technologie-Systeme GmbH.
 Aotal 28.
 Lassnitzhoehe A - 8301
 AUSTRIA

510(k) FIRM CONTACT: Al Lippincott
 Engineering Consulting Services, Inc.
 3150 E. 200th St.
 Prior Lake, MN 55372

TRADE NAME: FR.O.H. Calcaneus Repair System

COMMON NAME: Bone Plate & Cannulated Cancellous Bone Screw System

CLASSIFICATION: Plate, Fixation, Bone (see 21 CFR, Sec. 888.3030)
 Smooth or threaded metallic bone fixation fastener
 (see 21 CFR, Sec. 888.3040)

DEVICE PRODUCT CODE: HRS & HWC

SUBSTANTIALLY EQUIVALENT DEVICES Synthes Locking Calcaneal Plates (K991407)
 DePuy/Ace Calcaneal Peri-articular Plate (K993465)
 Ace/DePuy Cannulated Cancellous Bone Screw (K903810)

DEVICE DESCRIPTION: The I.T.S. FR.O.H. Calcaneus Repair System is a combination of fracture reduction and alignment instrumentation with either calcaneus plate and/or cannulated cancellous screw fixation across calcaneal heel bone fracture site(s). The 15-hole Calcaneus Plate (universal left and right in 2 sizes) is made from CP titanium according to ASTM F 67-00 and corresponding plate locking and self-tapping 3.5 & 4.2mm Cancellous Screws from 6-4 alloyed titanium according to ASTM F 136-02. The fully threaded 7.3mm Cannulated Cancellous Screw is in various lengths from 50 to 90mm in 5mm increment sizes and is also made from 6-4 alloyed titanium according to ASTM F 136-02. All titanium plates and screws are surface conditioned with a TIODIZE, Type II preparation.

201

510(k) Summary of Safety and Effectiveness Continued:


INTENDED USE: The I.T.S. FR.O.H. Calcaneus Repair System is used to stabilize a intra-articular and/or extra-articular fracture(s) and/or osteotomy of the calcaneus heel bone in the foot.

BASIS OF SUBSTANTIAL EQUIVALENCE: The I.T.S. FR.O.H. Calcaneus Repair System Calcaneal Plate is substantially equivalent to the Synthes (**K991407**) and, DePuy/Ace (**K993465**) plating systems.
The I.T.S. FR.O.H. Calcaneus Repair System Cannulated Cancellous Screw is substantially equivalent to the DePuy/Ace (**K903810**) cannulated screw system.

SUMMARY OF SAFETY AND EFFECTIVENESS: The I.T.S. FR.O.H. Calcaneus Repair System is shown to be safe and effective for use in fracture fixation of the calcaneus heel bone in the foot.



Implant-Technology-Systems



FR.O.H. Distraction

Calcaneus
Repair
Concept

Fixation by
angle stable
plate

according to


Univ. Prof. Dr. FRöhlich,
Prim. Dr. Ortner,
Prim. Univ. Prof. Dr. Hertz

Screwing

7.3 mm CANNULATED CANCELLOUS SCREW, fully threaded for fractures of Calcaneus

PROPERTIES

- Diameter Outside 7.3 mm
- Length from 50 to 90 mm in 5 mm Steps
- 3.2 mm Guide Wire applicable
- Selfcutting and selfdrilling
- Material: Ti Al6 V4
- Surface Treatment



IMPLANTS

- | | | |
|---|----------|--|
| 1 | 31731-xx | d=7.3 mm cannulated cancellous screw, fully threaded, Length=50-90mm, in 5mm Steps |
|---|----------|--|

INSTRUMENTS

- | | | |
|---|-----------|--|
| 2 | 35324-228 | d=3.2mm calibrated Guide Wire, width thread, L= 228 mm |
| 3 | 59321 | Hook Depth Gauge for calibrated Guide Wire |
| 4 | 53032 | Handle, flat |
| 5 | 54502-120 | SW 5.0mm cannulated Screw Driver Shank, L=120mm |
| 6 | 50172 | Sterilization Case |

CALCANEUS DISTRACTOR AND HEEL PRESS

INSTRUMENTS

- | | |
|-----------|-----------------------------------|
| 210530001 | Calcaneus Distractor, Fröhlich MD |
| 210530011 | Heel Press |
| 35324-228 | d=3.2mm Guide Wire, L= 228mm |



204

DEC 03 2001

K013044
page 1 of 1

3. **Summary of Safety and Effectiveness Information**

Sponsor	Synthes (USA) 1690 Russell Road Paoli, PA 19301
Company Contact	Matthew M. Hull (610) 647-9700 ext. 7191
Name of the Device	Synthes 3.9 mm Pelvic Screw
Device Classification(s)	Class II, §888.3040 – Screw, Fixation, Bone
Substantial Equivalence	Documentation was provided which demonstrated the Synthes 3.9 mm Pelvic Screw to be substantially equivalent to other legally marketed devices.
Device Description	The Synthes 3.9 mm Pelvic Screw is a machined metallic screw with a hex drive head.
Indications	The 3.9mm pelvic screw is indicated for fracture fixation of the pelvis as well as for periacetabular osteotomies.
Materials	Stainless Steel

MAR 20 2006

P. 1/2

K060156

Section XII: 510(k) Summary of Safety and Effectiveness

**SAFE MEDICAL DEVICES ACT OF 1990
510(k) Summary**

NAME OF FIRM: I.T.S. Implantat-Technologie-Systeme GmbH.
Aural 28.
Lassnitzhoehe A - 8301
AUSTRIA

510(k) FIRM CONTACT: Al Lippincott
Engineering Consulting Services, Inc.
3150 E. 200th St.
Prior Lake, MN 55372

TRADE NAME: Straight Plate with Angular Stability & Screw System

COMMON NAME: Bone Plate & Screw System

CLASSIFICATION: Plate, Fixation, Bone (see 21 CFR, Sec. 888.3030),
Screw, Fixation, Bone (see 21 CFR, Sec. 888.3040),
Washer, Bolt, Nut (see 21 CFR, Sec. 888.3030)

DEVICE PRODUCT CODE: HRS

SUBSEQUENT PRODUCT CODE: HWC, HTN

SUBSTANTIALLY I.T.S. GmbH Clavicularplate with Angular Stability (K050852)
Smith & Nephew Peri-Loc Locking Bone Plates and Locking Bone
Screw System (K051735)

EQUIVALENT DEVICES Zimmer Periarticular Locking Plates and Screws (K051098)
Synthes 4.5mm LCP Straight Reconstruction Plates (K051986)
Acumed Congruent Plate System (K012655)
I.T.S. GmbH FR.O.H. Calcaneus Repair System (K051643)
Synthes Sterile 3.5mm and 4.0mm Cannulated Screw (K963192)
Zimmer/Pioneer Cannulated Screw System (K003496)
Synthes (USA) Spherical Washers (K052483)

DEVICE DESCRIPTION: The I.T.S. Straight Plate with Angular Stability is a low-profile 4, 6, or 8 hole plate with various length cortical and/or cancellous self-tapping stabilization locking and/or compression screws. The Straight Plate is made from CP titanium according to ASTM F 67-00 and all screws are made from 6-4 Alloyed Titanium according to ASTM F 136-02. The plate and screws are surface conditioned with a TIODIZE, Type II preparation.

The I.T.S. Screw System is a group of cannulated fracture fixation screws in various diameters of 4.0mm, 6.5mm, and 7.3mm and lengths. A complement of flat and spherical Washers are available with the system. All screws and washers are made from 6-4 Alloyed Titanium according to ASTM F 136-02 and are surface conditioned with a TIODIZE, Type II preparation.

INTENDED USE:

The intended use of the I.T.S. Straight Plate with Angular Stability is to stabilize an osteotomy or fracture of small bones, long bones, the pelvis and the calcaneus in an adult or pediatric patient.

Indications for use include comminuted fractures, supercondylar fractures, intra-articular and extra-articular condylar fractures, fractures in osteopenic bone, nonunions, and malunions. And, as well, a fracture or osteotomy of the tibia, fibula, femoral condyle, acetabulum, humerus, ulna, middle hand and middle foot bones; treatment of the calcaneal; hip arthrodesis, and provisional hole fixation. The 3.5mm Cortical and Cancellous Angle Stable Screws used in conjunction with the Straight Plate may be used only on small bones.

The intended use of the I.T.S. Screw System is for corrective osteotomy or internal fracture fixation of the patella, pelvis, ankle, and long bones in an adult or pediatric patient.

For the 4.0mm Cannulated Cancellous Screw, indications for use are for radial and ulnar fractures, fractures of the proximal/distal humerus and of the patella, and for tendon fixation, mainonerve injuries and disruption of the syndesmosis with bimalleolar or supramalleolar fractures and the instability of the talus centering.

For the 6.5mm Cannulated Cancellous Screw, indications for use are for fractures of the femoral neck, tibiaplateau, of the sacrum and the articular cavity of the hip joint, and the metaphyseal fractures of the distal femur and distal tibia, fixation of the Ileo-sacral joint, and fusion of the foot and ank

For the 7.3mm Cannulated Cancellous Screw, indications for use are for fractures of the calcaneus, femoral neck, tibiaplateau, and of the sacrum and the articular cavity of the hip joint, fusion of the foot and ankle, fixation of the Ileo-sacral joint, and metaphyseal fractures of the distal femur and distal tibia.

The system(s) is not intended for spinal use.


BASIS OF SUBSTANTIAL EQUIVALENCE:

The I.T.S. Straight Plate with Angular Stability is substantially equivalent to the Smith & Nephew, Zimmer, Synthes, Acumed, and I.T.S. GmbH stabilizing bone plate systems. The I.T.S. Screw System is substantially equivalent to the I.T.S. GmbH, Zimmer, and Synthes cannulated screw and washer systems.

SUMMARY OF SAFETY AND EFFECTIVENESS:

The I.T.S. Straight Plate with Angular Stability and Screw System is shown to be safe and effective for use in fracture fixation of small and long bones in the body.

**Ø 7.3 MM CANNULATED CANCELLOUS-TENSION SCREW
WITH VARIABLE THREAD LENGTH**



Length (mm)	Item No.
40	31734-40
50	31734-50
55	31734-55
60	31734-60
65	31734-65
70	31734-70
75	31734-75
80	31734-80
85	31734-85
90	31734-90
95	31734-95
100	31734-100
105	31734-105
110	31734-110
115	31734-115
120	31734-120

PROPERTIES

- Back-tapping flank
- Ø 7.3 mm
- Length from 40 to 120 mm in 5 mm steps
- Ø 3.2 mm guide wire
- Self-drilling & self-tapping
- Material: Ti 6Al 4V ELI
- I.T.S. surface treatment



end-to-end thread



Length (mm)	Item No.
Screw	
40	31734-40
50	31734-50
55	31734-55
60	31734-60
65	31734-65
70	31734-70
75	31734-75
80	31734-80
85	31734-85
90	31734-90
95	31734-95
100	31734-100
105	31734-105
110	31734-110
115	31734-115
120	31734-120



Ø 3.2 mm calibrated guide wire, with thread L=225 mm
Depth gauge for calibrated guide wire
Hand grip (H)
SW Ø 9 mm cannulated allen head, L=120 mm

**Ø 7.3 MM CANNULATED CANCELLOUS-TENSION SCREW
WITH 15 MM THREAD LENGTH**



Length (mm)	Item No.
40	31733-40
50	31733-50
55	31733-55
60	31733-60
65	31733-65
70	31733-70
75	31733-75
80	31733-80
85	31733-85
90	31733-90
95	31733-95
100	31733-100
105	31733-105
110	31733-110
115	31733-115
120	31733-120

- Back-tapping flank
- Ø 7.3 mm
- Length from 40 to 120 mm in 5 mm steps
- Ø 3.2 mm guide wire
- Self-drilling & self-tapping
- Material: Ti 6Al 4V ELI
- I.T.S. surface treatment

Metatarsal fractures of the distal femur and distal tibia
Fractures of the acetabulum and the articular cavity of the hip joint
Fracture of the femoral neck
Fractures of the tibial plateau
Fracture of foot and ankle
Fracture of the scapula
Fracture of the femoral neck



**Ø 7.3 MM CANNULATED CANCELLOUS- TENSION SCREW
WITH END TO END THREAD**

INDICATIONS

- Fractures of the calcaneus
- Fractures of the femoral neck


On request also length from 120 mm in 5 mm steps available.

Length (mm)	Item No.
Screw	
50	31731-50
55	31731-55
60	31731-60
65	31731-65
70	31731-70
75	31731-75
80	31731-80
85	31731-85
90	31731-90
95	31731-95
100	31731-100
105	31731-105
110	31731-110
115	31731-115
120	31731-120

PROPERTIES

- Back-tapping flank
- Ø 7.3 mm
- length from 50 to 90 mm in 5 mm steps
- Ø 3.2 mm guide wire
- Self-drilling & self-tapping
- Material: Ti 6Al 4V ELI
- I.T.S. surface treatment

end-to-end thread



INSTRUMENTE

35324-228

Ø 3.2 mm calibrated guide wire, with thread, l= 228 mm

59321

Depth gauge for calibrated guide wire

53032

Hand grip, flat


54502-120

SW Ø 5.0 mm cannulated allen insert, l=120 mm

**Ø 7.3 MM CANNULATED CANCELLOUS- TENSION SCREW
WITH 16 MM THREAD LENGHT**

PROPERTIES

- Back-tapping flank
- Ø 7.3 mm
- Length from 80 to 120 mm in 5 mm steps
- Ø 3.2 mm guide wire
- Self-drilling & self-tapping
- Material: Ti 6Al 4V ELI
- I.T.S. surface treatment



Length (mm)	Item No.
Screw	
80	31732-80
85	31732-85
90	31732-90
95	31732-95
100	31732-100
105	31732-105
110	31732-110
115	31732-115
120	31732-120

INDICATIONS

- Fractures of the femoral neck
- Fractures of the tibiaplateau
- Fusion of foot and ankle
- Fixation of ileo-sacral joint
- Fractures of the sacrum and the articular cavity of the hip joint
- Metaphysal fractures of the distal femur and distal tibia

209



6.5 mm Cannulated Cancellous Screw

Thread 22 mm



PROPERTIES

- back-tapping flank
- 6.5 mm
- length from 25 to 120 mm in 5 mm steps
- 3.2 mm guide wire

INDICATIONS

- fractures of the femoral neck
- metaphysal fractures of the distal femur and distal tibia
- fractures of the tibial plateau
- fractures of the sacrum and the articular cavity of the hip joint
- fixation of the lileo-sacral joint
- fusion of foot and ankle


6.5 mm Cannulated SCFE Screw

Transcutaneous Screwing of the Slipped Capital Femoral Epiphysis



PROPERTIES

- 10 mm thread
- 3.2 mm guide wire
- Back-tapping flank
- Large countersink head for easy removableness (trocar 7.5 mm / 30°)



INDICATIONS

- loosening of the epiphysis femoralis

GAUGE 2.5
WITH CLAMP

7.3 mm Cannulated Cancellous Screw

fully threaded




PROPERTIES

- 7.3 mm
- length from 50 to 120 mm in 5 mm steps
- 3.2 mm guide wire

INDICATIONS

- fractures of the calcaneus
- fractures of the femoral neck



variable thread



PROPERTIES

- back-tapping flank
- 7.3 mm
- length from 40 to 120 mm in 5 mm steps
- 3.2 mm guide wire

INDICATIONS

- fractures of the femoral neck
- fractures of the tibial plateau
- fusion of foot and ankle
- fixation of lileo-sacral joint
- fractures of the sacrum and the articular cavity of the hip joint
- metaphysal fractures of the distal femur and distal tibia

16 mm thread



PROPERTIES

- back-tapping flank
- 7.3 mm
- length from 80 to 120 mm in 5 mm steps
- 3.2 mm guide wire

4.0 mm Cannulated Cancellous Screw

fully threaded




PROPERTIES

- 4.0 mm
- length from 45 to 60 mm in 5 mm steps
- 1.6 mm guide wire

INDICATIONS

- radial fractures
- ulnar fractures
- fractures of the proximal and distal humerus
- fractures of the patella
- tendon fixation

variable thread



PROPERTIES

- 4.0 mm
- length from 16 to 50 mm in 2 mm steps and 50 to 70 mm in 5 mm steps
- 1.6 mm guide wire

9 mm thread

designed for Dens Axis



PROPERTIES

- 4.0 mm
- length from 16 to 50 mm in 2 mm steps and 50 to 70 mm in 5 mm steps
- 1.6 mm guide wire

INDICATIONS

- fractures of the Dens Axis

4.0 MM WRENCH
TAPPING PIECE




4.0 MM GAUGE



Washers

Washer curved



Washer flat



for 7.3 mm / 4.0 mm
cancellous screw variable thread and
7.3 mm / 4.0 mm
cancellous screw 16 mm / 9 mm thread

FEATURES

- selftapping and selfdrilling
- material: TiAL6V4 ELI
- I.T.S. surface treatment

ALBERT LIPPINCOTT

From: "CDRH OC RLMAIL" <RLMAIL@CDRH.FDA.GOV>
To: "allippincott@msn.com" <allippincott@msn.com>
Sent: Friday, July 21, 2006 12:53 PM
Subject: Annual Establishment Registration for 2007 for 3004369035

Dear Official Correspondent:

This letter confirms that the following establishment is now registered for the entire calendar year of 2007:

I.T.S. IMPLANTAT-TECHNOLOGIE-SYSTEME GMBH
AUTAL 28
LASSNITZHOEHE,
AUSTRIA A-8301

The establishment listed above is registered as a medical device establishment, performing the following operation(s):

MANUFACTURER

Please be aware that this confirmation does not indicate or imply approval of any product or activity at the establishment. It is also not a license or certification of any kind.

Please remember to report any changes to your establishment registration information using the form FDA 2891, Registration of Device Establishment, within the timeframes specified in Title 21 of the Code of Federal Regulations, Part 807.

If you have any questions regarding your registration or device listings, please send an email to reglist@cdrh.fda.gov or phone 240-276-0111.

Food and Drug Administration
Center for Devices and Radiological Health, HFZ-308
9200 Corporate Blvd.
Rockville, MD 20850-4015

212

7/21/2006

DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service
Food and Drug Administration
Memorandum

From: Reviewer(s) - Name(s) N. K. MISHRA

Subject: 510(k) Number K063/66

To: The Record - It is my recommendation that the subject 510(k) Notification:

- Refused to accept.
- Requires additional information (other than refuse to accept).
- Is substantially equivalent to marketed devices.
- NOT substantially equivalent to marketed devices.
- Other (e.g., exempt by regulation, not a device, duplicate, etc.)

- Is this device subject to Section 522 Postmarket Surveillance? YES NO
- Is this device subject to the Tracking Regulation? YES NO
- Was clinical data necessary to support the review of this 510(k)? YES NO
- Is this a prescription device? YES NO
- Was this 510(k) reviewed by a Third Party? YES NO
- Special 510(k)? YES NO
- Abbreviated 510(k)? Please fill out form on H Drive 510k/boilers YES NO

- Truthful and Accurate Statement Requested Enclosed
- A 510(k) summary OR A 510(k) statement
- The required certification and summary for class III devices
- The indication for use form

Combination Product Category (Please see algorithm on H drive 510k/Boilers) N

Animal Tissue Source YES NO Material of Biological Origin YES NO

The submitter requests under 21 CFR 807.95 (doesn't apply for SEs):

- No Confidentiality
- Confidentiality for 90 days
- Continued Confidentiality exceeding 90 day:

Predicate Product Code with class: Additional Product Code(s) with panel (optional):

HRS/OR/II 21CFR§ 888.3030 JDA/OR/II 21CFR 888.3030

Review: [Signature] JDA HWC/OR/II 21CFR 888.3040
(Branch Chief) (Branch Code) (Date) 12/22/06

Final Review: [Signature] 12/22/06
(Division Director) (Date)

Internal Administrative Form

	YES	NO
1. Did the firm request expedited review?		X
2. Did we grant expedited review?		X
3. Have you verified that the Document is labeled Class III for GMP purposes?		X
4. If, not, has POS been notified?		
5. Is the product a device?	X	
6. Is the device exempt from 510(k) by regulation or policy?	X	X
7. Is the device subject to review by CDRH?		
8. Are you aware that this device has been the subject of a previous NSE decision?		Y
9. If yes, does this new 510(k) address the NSE issue(s), (e.g., performance data)?		
10. Are you aware of the submitter being the subject of an integrity investigation?		X
11. If, yes, consult the ODE Integrity Officer.		
12. Has the ODE Integrity Officer given permission to proceed with the review? (Blue Book Memo #I91-2 and Federal Register 90N0332, September 10, 1991.		

**SCREENING CHECKLIST
FOR ALL PREMARKET NOTIFICATION [510(k)] SUBMISSIONS**

510(k) Number 063166

The cover letter clearly identifies the type of 510(k) submission as (Check the appropriate box):

- Special 510(k) - Do Sections 1 and 2
- Abbreviated 510(k) - Do Sections 1, 3 and 4
- Traditional 510(k) or no identification provided - Do Sections 1 and 4

Section 1: Required Elements for All Types of 510(k) submissions:

	Present or Adequate	Missing or Inadequate
Cover letter, containing the elements listed on page 3-2 of the Premarket Notification [510] Manual.	X	
Table of Contents.	X	
Truthful and Accurate Statement.	X	
Device's Trade Name, Device's Classification Name and Establishment Registration Number.	X	
Device Classification Regulation Number and Regulatory Status (Class I, Class II, Class III or Unclassified).	Y	
Proposed Labeling including the material listed on page 3-4 of the Premarket Notification [510] Manual.	X	
Statement of Indications for Use that is on a separate page in the premarket submission.	X	
Substantial Equivalence Comparison, including comparisons of the new device with the predicate.	X	
510(k) Summary or 510(k) Statement.	X	
Description of the device (or modification of the device) including diagrams, engineering drawings, photographs or service manuals.	X	
Identification of legally marketed predicate device. *	X	
Compliance with performance standards. * [See Section 514 of the Act and 21 CFR 807.87 (d).]	X	
Class III Certification and Summary. **		
Financial Certification or Disclosure Statement for 510(k) notifications with a clinical study. * [See 21 CFR 807.87 (i)]		
510(k) Kit Certification ***		

* - May not be applicable for Special 510(k)s.

** - Required for Class III devices, only.

*** - See pages 3-12 and 3-13 in the Premarket Notification [510] Manual and the Convenience Kits Interim Regulatory Guidance.

Section 2: Required Elements for a SPECIAL 510(k) submission:

	Present	Inadequate or Missing
Name and 510(k) number of the submitter's own, unmodified predicate device.		
A description of the modified device and a comparison to the sponsor's predicate device.		
A statement that the intended use(s) and indications of the modified device, as described in its labeling are the same as the intended uses and indications for the submitter's unmodified predicate device.		
Reviewer's confirmation that the modification has not altered the fundamental scientific technology of the submitter's predicate device.		
A Design Control Activities Summary that includes the following elements (a-c):		
a. Identification of Risk Analysis method(s) used to assess the impact of the modification on the device and its components, and the results of the analysis.		
b. Based on the Risk Analysis, an identification of the required verification and validation activities, including the methods or tests used and the acceptance criteria to be applied.		
c. A Declaration of Conformity with design controls that includes the following statements:		
A statement that, as required by the risk analysis, all verification and validation activities were performed by the designated individual(s) and the results of the activities demonstrated that the predetermined acceptance criteria were met. This statement is signed by the individual responsible for those particular activities.		
A statement that the manufacturing facility is in conformance with the design control procedure requirements as specified in 21 CFR 820.30 and the records are available for review. This statement is signed by the individual responsible for those particular activities.		

Section 3: Required Elements for an ABBREVIATED 510(k)* submission:

	Present	Inadequate or Missing
For a submission, which relies on a guidance document and/or special control(s), a summary report that describes how the guidance and/or special control(s) was used to address the risks associated with the particular device type. (If a manufacturer elects to use an alternate approach to address a particular risk, sufficient detail should be provided to justify that approach.)		
For a submission, which relies on a recognized standard, a declaration of conformity [For a listing of the required elements of a declaration of conformity, SEE Required Elements for a Declaration of Conformity to a Recognized Standard, which is posted with the 510(k) boilers on the H drive.]		

For a submission, which relies on a recognized standard without a declaration of conformity, a statement that the manufacturer intends to conform to a recognized standard and that supporting data will be available before marketing the device.		
For a submission, which relies on a non-recognized standard that has been historically accepted by FDA, a statement that the manufacturer intends to conform to a recognized standard and that supporting data will be available before marketing the device.		
For a submission, which relies on a non-recognized standard that has <u>not</u> been historically accepted by FDA, a statement that the manufacturer intends to conform to a recognized standard and that supporting data will be available before marketing the device <u>and</u> any additional information requested by the reviewer in order to determine substantial equivalence.		
Any additional information, which is not covered by the guidance document, special control, recognized standard and/or non-recognized standard, in order to determine substantial equivalence.		

- * - When completing the review of an abbreviated 510(k), please fill out an Abbreviated Standards Data Form (located on the H drive) and list all the guidance documents, special controls, recognized standards and/or non-recognized standards, which were noted by the sponsor.

Section 4: Additional Requirements for ABBREVIATED and TRADITIONAL 510(k) submissions (If Applicable):

	Present	Inadequate or Missing
a) Biocompatibility data for all patient-contacting materials, OR certification of identical material/formulation:		
b) Sterilization and expiration dating information:	X	
i) sterilization process		
ii) validation method of sterilization process		
iii) SAL		
iv) packaging		
v) specify pyrogen free		
vi) ETO residues		
vii) radiation dose		
viii) Traditional Method or Non-Traditional Method		
c) Software Documentation:		

Items with checks in the "Present or Adequate" column do not require e additional information from the sponsor. Items with checks in the "Missing or Inadequate" column must be submitted before substantive review of the document.

Passed Screening Yes No
 Reviewer: Wishu 12/15/06
 Concurrence by Review Branch: Jonita
 Date: 12/20/06

(b) (4)



510(k) Submission
I.T.S. Implantat-Technologie-Systeme GmbH
Pelvic Reconstruction System (PRS)

Section X: Labeling

See attached sample labeling for various components of the I.T.S. Pelvic Reconstruction System (PRS) designating the manufacturer and address, product catalog number, product description, quantity in package, material, lot/production number traceability, and sterilization designation as **NONSTERILE**.

A package insert for the I.T.S. Pelvic Reconstruction System (PRS) is currently in development that will disclose the Indications, Warnings, Precautions, Contraindications, Sterilization as NONSTERILE as labeled (with recommendations for in-hospital gravity autoclave per AAMI guidelines), and Materials. A clause will be included stating: "The system(s) has not been studied in spinal use, and is not intended for use in vertebral column fracture or fusion procedures. Also, the statement "Federal Law (USA) restricts this device to sale by or on order of a board certified physician" will be shown.

3CB 12/22/06

THE 510(K) DOCUMENTATION FORMS ARE AVAILABLE ON THE LAN UNDER 510(K) BOILERPLATES TITLED "DOCUMENTATION" AND MUST BE FILLED OUT WITH EVERY FINAL DECISION (SE, NSE, NOT A DEVICE, ETC.).

"SUBSTANTIAL EQUIVALENCE" (SE) DECISION MAKING DOCUMENTATION

K 063166

Reviewer: N.K. MISHRA

Division/Branch: DGRND

Device Name: ITS Pelvic Reconstruction SET

Product To Which Compared (510(K) Number If Known): K0017270 K031573
K001614

	YES	NO	
1. Is Product A Device	<input checked="" type="checkbox"/>	<input type="checkbox"/>	If NO = Stop
2. Is Device Subject To 510(k)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	If NO = Stop
3. Same Indication Statement?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	If YES = Go To 5
4. Do Differences Alter The Effect Or Raise New Issues of Safety Or Effectiveness?	<input type="checkbox"/>	<input type="checkbox"/>	If YES = Stop NE
5. Same Technological Characteristics?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	If YES = Go To 7
6. Could The New Characteristics Affect Safety Or Effectiveness?	<input type="checkbox"/>	<input type="checkbox"/>	If YES = Go To 8
7. Descriptive Characteristics Precise Enough?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	If NO = Go To 10 If YES = Stop SE
8. New Types Of Safety Or Effectiveness Questions?	<input type="checkbox"/>	<input type="checkbox"/>	If YES = Stop NE
9. Accepted Scientific Methods Exist?	<input type="checkbox"/>	<input type="checkbox"/>	If NO = Stop NE
10. Performance Data Available?	<input type="checkbox"/>	<input type="checkbox"/>	If NO = Request Data
11. Data Demonstrate Equivalence?	<input type="checkbox"/>	<input type="checkbox"/>	Final Decision:

Note: In addition to completing the form on the LAN, "yes" responses to questions 4, 6, 8, and 11, and every "no" response requires an explanation.


1. Intended Use:
2. Device Description: Provide a statement of how the device is either similar to and/or different from other marketed devices, plus data (if necessary) to support the statement. Is the device life-supporting or life sustaining? Is the device implanted (short-term or long-term)? Does the device design use software? Is the device sterile? Is the device for single use? Is the device over-the-counter or prescription use? Does the device contain drug or biological product as a component? Is this device a kit? Provide a summary about the devices design, materials, physical properties and toxicology profile if important.

EXPLANATIONS TO "YES" AND "NO" ANSWERS TO QUESTIONS ON PAGE 1 AS NEEDED

1. Explain why not a device:
2. Explain why not subject to 510(k):
3. How does the new indication differ from the predicate device's indication:
4. Explain why there is or is not a new effect or safety or effectiveness issue:
5. Describe the new technological characteristics:
6. Explain how new characteristics could or could not affect safety or effectiveness:
7. Explain how descriptive characteristics are not precise enough:
8. Explain new types of safety or effectiveness questions raised or why the questions are not new:
9. Explain why existing scientific methods can not be used:
10. Explain what performance data is needed:
11. Explain how the performance data demonstrates that the device is or is not substantially equivalent:

ATTACH ADDITIONAL SUPPORTING INFORMATION

510(k) "SUBSTANTIAL EQUIVALENCE" DECISION-MAKING PROCESS



- * 510(k) Submissions compare new devices to marketed devices. FDA requests additional information if the relationship between marketed and "predicate" (pre-Amendments or reclassified post-Amendments) devices is unclear.
- ** This decision is normally based on descriptive information alone, but limited testing information is sometimes required.
- *** Data maybe in the 510(k), other 510(k)s, the Center's classification files, or the literature.

5