

K122753



DEC 14 2012

## 510(k) Summary

Page 5-1

**Preparation Date: August 28, 2012**

### Company Information:

American Orthodontics  
 1714 Cambridge Avenue  
 Sheboygan, WI 53081  
 Phone: 920-457-5051  
 Fax: 920-457-5773

### Contact Information:

Trang Adams / Regulatory Affairs Specialist  
 1536 N. 18<sup>th</sup> Street  
 Sheboygan, WI 53081  
 Phone: 920-457-5051 Ext. 4251  
 Fax: 920-457-5773  
 E-Mail: tadams@americanortho.com

### Device Information:

Trade Name: Empower Clear  
 Common Name: Ceramic Brackets  
 Classification Name: Bracket, Ceramic, Orthodontic  
 Classification Code: NJM  
 Regulation Number: 872.5470

### Equivalent Legally Marketed Devices Information:

<u>510(k) #</u>	<u>Product Name</u>	<u>Device Manufacturer</u>
K080749	Radiance	American Orthodontics
K060837	In-Ovation C	Densply

### Description of the Device:

The Empower Clear line of products is single-use devices intended for use in conjunction with comprehensive orthodontics to control the movement of individual teeth. The Ceramic Bracket combines the aesthetics of a ceramic bracket with the versatility and ease of self ligation.

These brackets are comprised of several geometries that vary from bracket to bracket, corresponding to the intended tooth. These geometries contribute to the fit of the bracket to the tooth and also impart the axial control of the energy from the archwire.

### Indications for Use:

Empower Clear brackets are intended for orthodontic movement of teeth as diagnosed by an orthodontist. It is used temporarily and is removed upon completion of orthodontic treatment. Empower Clear brackets are intended to be single use only.

**Technological Characteristics Information:**

The material, Alumina Oxide [Al<sub>2</sub>O<sub>3</sub>], is used in the manufacturing of the Radiance brackets and In-Ovation C brackets – which is the same material used for Empower Clear.

The In-Ovation C clip material is the same as the Empower Clear clip material. The Radiance bracket does not have a clip.

The function and performance of the Empower Clear brackets are substantially equivalent to the predicate devices, as outlined in the following table.

Product Parameter	Device Name / Manufacturer			Substantial Equivalence Analysis
	Radiance / American Orthodontics	In-Ovation / Dentsply	Empower / American Orthodontics	
510(k) Number	K080749	K060837	Pending	N/A
Material	Al <sub>2</sub> O <sub>3</sub>	Al <sub>2</sub> O <sub>3</sub>	Al <sub>2</sub> O <sub>3</sub>	Equivalent
Intended Use	Orthodontic treatment is used to correct dental deficiencies and to improve the appearance of the patient. The brackets, arch wire and elastic o-rings form a force system that is designed to gradually move teeth into a normal alignment.	The Innovation C is intended for orthodontic movement of natural teeth, excluding the mandibular bicuspid teeth.	Ceramic Brackets are intended for orthodontic movement of teeth as diagnosed by an orthodontist. It is used temporarily and is removed upon completion of orthodontic treatment. Ceramic Brackets are intended to be single use only.	Equivalent
Single Use	YES	YES	YES	Equivalent
Non-Sterile Packaging	YES	YES	YES	Equivalent

**Biocompatibility Testing:**

Biocompatibility testing conducted on Empower Clear brackets indicates that the bracket material is safe for use. The bracket material, Aluminum Oxide [Al<sub>2</sub>O<sub>3</sub>], was found to be free of harmful extractables. No oral mucosa irritation or skin sensitization was detected with the material.

**Summary:**

The function and performance of Empower Clear bracket is similar to the predicates. There are no changes in the intended use and fundamental scientific technology. All of the materials used in the device have been used in legally marketed American Orthodontics devices. Minor differences in



technological characteristics do not raise new types of safety and effectiveness questions.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

December 14, 2012

Mr. Trang Adams  
Regulatory Affairs Specialist  
American Orthodontics  
1536 North 18<sup>th</sup> Street  
SHEBOYGAN WI 53081

Re: K122753  
Trade/Device Name: Empower Clear  
Regulation Number: 21 CFR 872.5470  
Regulation Name: Orthodontic Plastic Bracket  
Regulatory Class: II  
Product Code: NJM  
Dated: October 23, 2012  
Received: October 26, 2012

Dear Mr. Adams:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Page 2 – Mr. Adams

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Susan Runner DDS, MA

2012.12.14

10:46:01 -05'00'

Anthony D. Watson, B.S., M.S., M.B.A.  
Director  
Division of Anesthesiology, General Hospital,  
Respiratory, Infection Control and  
Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

K 122753



AMERICAN  
ORTHODONTICS

### Indications for Use Statement

510(k) Number (if known): K 122753  
Unknown

Device Name:  
Empower Clear

**Indications for Use:**

Empower Clear brackets are intended for orthodontic movement of teeth as diagnosed by an orthodontist. It is used temporarily and is removed upon completion of orthodontic treatment. Empower Clear brackets are intended to be single use only.

**Prescription Use And/Or Over-The-Counter Use:**

- Prescription use by orthodontist only
- Not available Over-The-Counter [OTC]

2012.12.1

Susan Runner DDS, MA 4 10:47:20

-05'00'

**(Division Sign-Off)**  
**Division of Anesthesiology, General Hospital**  
**Infection Control, Dental Devices**

510(k) Number: K122753

K122753

v.1



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

December 14, 2012

Mr. Trang Adams  
Regulatory Affairs Specialist  
American Orthodontics  
1536 North 18<sup>th</sup> Street  
SHEBOYGAN WI 53081

Re: K122753  
Trade/Device Name: Empower Clear  
Regulation Number: 21 CFR 872.5470  
Regulation Name: Orthodontic Plastic Bracket  
Regulatory Class: II  
Product Code: NJM  
Dated: October 23, 2012  
Received: October 26, 2012

Dear Mr. Adams:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Page 2 – Mr. Adams

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucml15809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Susan Runner DDS, MA

2012.12.14

10:46:01 -05'00'

Anthony D. Watson, B.S., M.S., M.B.A.  
Director  
Division of Anesthesiology, General Hospital,  
Respiratory, Infection Control and  
Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure



Page 3 – Mr. Adams

**Concurrence & Template History Page**  
 [THIS PAGE IS INCLUDED IN IMAGE COPY ONLY]



Full Submission Number: K122753

For Office of Compliance Contact Information:

[http://insideportlets.fda.gov:9010/portal/page?\\_pageid=197,415881&\\_dad=portal&\\_schema=PORTAL&org=318](http://insideportlets.fda.gov:9010/portal/page?_pageid=197,415881&_dad=portal&_schema=PORTAL&org=318)

For Office of Surveillance and Biometrics Contact Information:

[http://insideportlets.fda.gov:9010/portal/page?\\_pageid=197,415881&\\_dad=portal&\\_schema=PORTAL&org=423](http://insideportlets.fda.gov:9010/portal/page?_pageid=197,415881&_dad=portal&_schema=PORTAL&org=423)

<b>Digital Signature Concurrence Table</b>	
Reviewer Sign-Off	
Branch Chief Sign-Off	2012.12.14 Susan Runner DDS, MA 10:46:21 -05'00' 
Division Sign-Off	2012.12.14 Susan Runner DDS, MA 12:13:21 -05'00' 

Template Name: K1(A) – SE after 1996

Template History:

Date of Update	By	Description of Update
7/27/09	Brandi Stuart	Added Updates to Boiler Table
8/7/09	Brandi Stuart	Updated HFZ Table
1/11/10	Diane Garcia	Liability/Warranty sentence added at bottom of 1 <sup>st</sup> page
10/4/11	M. McCabe Janicki	Removed IFU sheet and placed in Forms
9/25/12	Edwena Jones	Added digital signature format
12/12/12	M. McCabe Janicki	Added an extra line between letter signature block and the word "Enclosure". Also, added a missing digit in 4-digit extension on letterhead zip code: "002" should be "0002".

K 122753



AMERICAN  
ORTHODONTICS

### Indications for Use Statement

510(k) Number (if known): K 122753  
Unknown

Device Name:  
Empower Clear

Indications for Use:  
Empower Clear brackets are intended for orthodontic movement of teeth as diagnosed by an orthodontist. It is used temporarily and is removed upon completion of orthodontic treatment. Empower Clear brackets are intended to be single use only.

**Prescription Use And/Or Over-The-Counter Use:**

- Prescription use by orthodontist only
- Not available Over-The-Counter [OTC]

2012.12.1

Susan Runner DDS, MA 4 10:47:20

-05'00'

**(Division Sign-Off)**  
Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K122753

\* \* \* COMMUNICATION RESULT REPORT ( DEC. 14. 2012 2:39PM ) \* \* \*

FAX HEADER 1:  
FAX HEADER 2:TRANSMITTED/STORED : DEC. 14. 2012 2:34PM  
MODE OPTION

ADDRESS

RESULT

PAGE

1887 MEMORY TX

919204575773

OK

3/3

REASON FOR ERROR  
E-1) HANG UP OR LINE FAIL  
E-3) NO ANSWERE-2) BUSY  
E-4) NO FACSIMILE CONNECTION

DEPARTMENT OF HEALTH &amp; HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

December 14, 2012

Mr. Trang Adams  
Regulatory Affairs Specialist  
American Orthodontics  
1536 North 18<sup>th</sup> Street  
SHEBOYGAN WI 53081

Re: K122753  
Trade/Device Name: Empower Clear  
Regulation Number: 21 CFR 872.5470  
Regulation Name: Orthodontic Plastic Bracket  
Regulatory Class: II  
Product Code: NJM  
Dated: October 23, 2012  
Received: October 26, 2012

Dear Mr. Adams:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.



**COVER SHEET MEMORANDUM**

**From:** Reviewer Name Mary E Brown  
**Subject:** 510(k) Number K122753/S001  
**To:** The Record

Please list CTS decision code SE

- Refused to accept (Note: this is considered the first review cycle, See Screening Checklist [http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0\\_5631/Screening%20Checklist%207%202%2007.doc](http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0_5631/Screening%20Checklist%207%202%2007.doc))
- Hold (Additional Information or Telephone Hold).
- Final Decision (SE, SE with Limitations, NSE (select code below), Withdrawn, etc.).

Not Substantially Equivalent (NSE) Codes

- NO NSE for lack of predicate
- NI NSE for new intended use
- NQ NSE for new technology that raises new questions of safety and effectiveness
- NU NSE for new intended use AND new technology raising new questions of safety and effectiveness
- NP NSE for lack of performance data
- NS NSE no response
- NL NSE for lack of performance data AND no response
- NM NSE pre-amendment device call for PMAs (515i)
- NC NSE post-amendment device requires PMAs
- NH NSE for new molecular entity requires PMA
- TR NSE for transitional device






Please complete the following for a final clearance decision (i.e., SE, SE with Limitations, etc.):		YES	NO
Indications for Use Page	<i>Attach IFU</i>	✓	
510(k) Summary /510(k) Statement	<i>Attach Summary</i>	✓	
Truthful and Accurate Statement.	<i>Must be present for a Final Decision</i>	✓	
Is the device Class III?		✓	
If yes, does firm include Class III Summary?	<i>Must be present for a Final Decision</i>		✓
Does firm reference standards? (If yes, please attach form from <a href="http://www.fda.gov/opacom/morechoices/fdaforms/FDA-3654.pdf">http://www.fda.gov/opacom/morechoices/fdaforms/FDA-3654.pdf</a> )		✓	
Is this a combination product? (Please specify category _____, see <a href="http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0_413b/CO-MBINATION%20PRODUCT%20ALGORITHM%20(REVISED%203-12-03).DOC">http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0_413b/CO-MBINATION%20PRODUCT%20ALGORITHM%20(REVISED%203-12-03).DOC</a> )			✓
Is this a reprocessed single use device? (Guidance for Industry and FDA Staff – MDUFMA - Validation Data in 510(k)s for Reprocessed Single-Use Medical Devices, <a href="http://www.fda.gov/cdrh/ode/guidance/1216.html">http://www.fda.gov/cdrh/ode/guidance/1216.html</a> )			✓
Is this device intended for pediatric use only?			✓
Is this a prescription device? (If both prescription & OTC, check both boxes.)		✓	
Did the application include a completed FORM FDA 3674, Certification with Requirements of ClinicalTrials.gov Data Bank?			✓
Is clinical data necessary to support the review of this 510(k)?			✓
For United States-based clinical studies <b>only</b> : Did the application include a completed FORM FDA 3674, Certification with Requirements of ClinicalTrials.gov Data Bank? (If study was			✓

conducted in the United States, and FORM FDA 3674 was not included or incomplete, then applicant must be contacted to obtain completed form.)	✓
Does this device include an Animal Tissue Source?	✓
All Pediatric Patients age ≤ 21	✓
Neonate/Newborn (Birth to 28 days)	✓
Infant (29 days - < 2 years old)	✓
Child (2 years - < 12 years old)	✓
Adolescent (12 years - < 18 years old)	✓
Transitional Adolescent A (18 - < 21 years old) Special considerations are being given to this group, different from adults age ≥ 21 (different device design or testing, different protocol procedures, etc.)	✓
Transitional Adolescent B (18 - ≤ 21; No special considerations compared to adults ⇒ 21 years old)	✓
Nanotechnology	✓
Is this device subject to the Tracking Regulation? (Medical Device Tracking Guidance, <a href="http://www.fda.gov/cdrh/comp/guidance/169.html">http://www.fda.gov/cdrh/comp/guidance/169.html</a> )	✓
Contact OC.	


Regulation Number	Class*	Product Code
872-5470	C. II	76 NJM

(\*If unclassified, see 510(k) Staff)

Additional Product Codes:

Review:	 (Branch Chief)	 (Branch Code)	 (Date)
Final Review:	 (Division Director)		 (Date)

### 510(k) "SUBSTANTIAL EQUIVALENCE" DECISION-MAKING PROCESS



\* 510(k) Submissions compare new devices to marketed devices. FDA requests additional information if the relationship between marketed and "predicate" (pre-Amendments or reclassified post-Amendments) devices is unclear.

\*\* This decision is normally based on descriptive information alone, but limited testing information is sometimes required.

\*\*\* Data may be in the 510(k), other 510(k)s, the Center's classification files, or the literature.



**DEPARTMENT OF HEALTH AND HUMAN SERVICES**

**MEMORANDUM**

Food and Drug Administration  
Office of Device Evaluation  
9200 Corporate Boulevard  
Rockville, MD 20850

**Premarket Notification [510(k)] Review  
Traditional/Abbreviated**

**K122753**

**Date:** November 30, 2012  
**To:** The Record  
**From:** Myra E. Browne, M.S., Biologist  
**Office/Division:** ODE/DAGRID  
**510(k) Holder:** American Orthodontics  
**Device Name:** Empower Clear  
**Contact:** Ms. Trang Adams  
**Phone:** 920-457-5051  
**Fax:** 920-457-5773  
**Email:** [tadams@americanortho.com](mailto:tadams@americanortho.com)

**Purpose and Submission Summary**

The 510(k) holder would like to introduce Empower Clear into interstate commerce.

Empower Clear is a ceramic orthodontic bracket that is intended to reposition teeth during orthodontic treatment.

Empower Clear is substantially equivalent (SE) to legally marketed ceramic orthodontic brackets because the information submitted by American Orthodontics, demonstrates that the device has the same indication and technological characteristics as legally marketed ceramic orthodontic brackets.

**Administrative Requirements**

	Yes	No	N/A
Indications for Use page (Indicate if: Prescription or OTC)	X		
Truthful and Accuracy Statement	X		
510(k) Summary or 510(k) Statement	X		
Standards Form	X		

**Indications for Use**

Empower Clear ceramic orthodontic brackets are intended for attachment to the patient's

teeth during orthodontic treatment to reposition the teeth to achieve correct dental deficiencies.

The indication of Empower Clear does not differ from that of legally marketed orthodontic ceramic brackets.

Device Description/Formulation

	Yes	No	N/A
Is the device life-supporting or life sustaining?		X	
Is the device an implant (implanted longer than 30 days)?		X	
Does the device design use software?		X	
Is the device sterile?		X	
Is the device reusable (not reprocessed single use)?		X	
Are "cleaning" instructions included for the end user?		X	

The purpose of this 510(K) is to introduce a new product to market. No novel features have been introduced.

Empower Clear is an orthodontic bracket line of products that are used in conjunction with comprehensive orthodontics to control the movement of individual teeth. The ceramic bracket is a self ligating bracket system. The Empower Bracket line consists of Empower Self-ligating stainless steel brackets, Empower Reconvertible stainless steel buccal tubes and Empower Clear brackets. Empower Clear is available as interactive or passive self-ligating systems. The interactive self-ligating system offers less friction at the beginning of treatment and more control in the later stages. The passive self-ligating system offers less friction and lower force mechanics throughout treatment, no wires are actively engage by the clip. The Empower Clear brackets are made of aluminum oxide which gives the brackets a clear appearance which are desirable for esthetics.

The Empower Clear bracket system consists of 20 brackets: a group of 4 identical brackets (upper bicuspids), a second group of 4 identical brackets (lower centrals and laterals) and 12 unique brackets. Slot sizes are available in 0.018" and 0.022". The cuspid and bicuspid brackets are available with and without hooks that are used for auxiliary attachments. The three axis force system for each bracket is defined by the torque, rotation and angulation component values.

Empower Clear orthodontic brackets use a colored dot identification system to facilitate placement of the brackets on the teeth via the location of one or two colored dots on the bracket tie wings. Each colored dot is approximately 0.015" in diameter. These dots are made of food coloring purchased by (b) (4) and are brushed away at the patients' initial brushings after placement. In addition to these dots the brackets contain Visual Positioning Aids in arch wire slots to assist wit the bracket alignment during bonding. These Visual Positioning Aids are removed with a scaler prior to final placement of the brackets so that the wire can be placed.



The chemical composition of Empower Clear is as follows:

Chemical	%
Aluminum oxide	99.99
Impurities which may exist as: Si, Fe, Na, K, Ca and/or Mg	0.01

(b) (4)

#### Contact History

The reviewer contacted the submitter via telephone on October 19, 2012, and spoke to Ms. Trang Adams to request additional information regarding the chemical composition, 510(k) summary, physical properties, removal of the brackets, etc. The additional information was received on October 26, 2012. The company did not submit the revised 510(k) summary with the biocompatibility information included. The revised summary was received on December 12, 2012. There are no outstanding issues.

#### Deficiencies

No deficiencies have been identified.

#### Labeling

The labeling for Empower Clear has been provided which includes instructions for use and an appropriate prescription statement as required by CFR 21.801.109. No unsubstantiated claims are purported. Labeling includes the proper bonding and debonding instructions.

#### Sterilization/Shelf Life/Reuse:

Empower Clear will be provided non-sterile and is not intended to be sterilized before use.

#### Biocompatibility

The formulation of Empower Clear includes no new components. This basic formulation is known to be biocompatible for this intended use. The company did submit cytotoxicity testing, oral mucosa testing and sensitization testing. The results of these tests demonstrate that Empower Clear is acceptable for its intended use as a ceramic orthodontic bracket.

#### Software

Empower Clear contains no software.

Electromagnetic Compatibility and Electrical, Mechanical and Thermal Safety

Empower Clear is not a mechanical or electrical device. Therefore, mechanical safety, electrical safety, EMC, and thermal safety are not applicable.

Performance Testing - Bench

Engineering performance test results are provided in Section, Device Comparison.

Performance Testing - Animal

Animal test results were not provided for Empower Clear. This type of information is not needed for the assessment of safety and effectiveness of this product.

Performance Testing - Clinical

Human test results were not provided for Empower Clear. This type of information is not needed for the assessment of safety and effectiveness of this product.

Device Comparison

Predicate Device: Radiance Ceramic Orthodontic Bracket (K080749) of American Orthodontics.

Physical Property	(b) (4)
Tensile strength (MPa)	
Density (grams/cubic centimeter)	
Melting point	
Vickers Hardness	
Transparency	

Empower Clear orthodontic brackets are comparable to other legally marketed orthodontic ceramic brackets on the market, especially the Radiance Orthodontic Bracket (K080749) also manufactured by American Orthodontics. Both Radiance and Empower Clear brackets are made of 99.99% aluminum oxide, and both have (b) (4)

The difference between the two products is that the Radiance material is monocrystalline, whereas Empower Clear is polycrystalline. In addition, Radiance brackets do not have a clip, and Empower Clear brackets have a clip. Both brackets have the visual placements aids for positioning and alignment, the Radiance aid is colored + and the Empower's aid is only in the slot with a -.

The physical properties of Empower Clear appear to be adequate for its intended use.

No new technological characteristics have been introduced in Empower Clear that could affect its safety or effectiveness.

**Substantial Equivalence Discussion**

	Yes	No
1. Same Indication Statement?	X	If YES = Go To 3
2. Do Differences Alter The Effect Or Raise New Issues of Safety Or Effectiveness?		If YES = Stop NSE
3. Same Technological Characteristics?	X	If YES = Go To 5
4. Could The New Characteristics Affect Safety Or Effectiveness?		If YES = Go To 6
5. Descriptive Characteristics Precise Enough?	X	If NO = Go To 8 If YES = Stop SE
6. New Types Of Safety Or Effectiveness Questions?		If YES = Stop NSE
7. Accepted Scientific Methods Exist?		If NO = Stop NSE
8. Performance Data Available?		If NO = Request Data
9. Data Demonstrate Equivalence?		Final Decision: SE

**Recommendation**

Regulation Number: 21 CFR 872.5470  
 Regulation Name: Orthodontic Ceramic Bracket  
 Regulatory Class: Class II  
 Product Code: NJM

**Myra E.  
Browne**

Digitally signed by Myra E. Browne  
 DN: c=US, o=U.S. Government, ou=HHS,  
 ou=FDA, ou=People, cn=Myra E. Browne,  
 0.9.2342.19200300.100.1.1=1300013790  
 Date: 2012.12.14 10:30:45 -05'00'

Myra E. Browne, M. S., Biologist  
 Reviewer

\_\_\_\_\_  
 Date

Susan Runner DDS, MA 2012.12.14  
 10:46:55 -05'00'

M. Susan Runner, DDS  
 Branch Chief

\_\_\_\_\_  
 Date



## 510(k) Summary

Page 5-1

Preparation Date: August 28, 2012

### Company Information:

American Orthodontics  
1714 Cambridge Avenue  
Sheboygan, WI 53081  
Phone: 920-457-5051  
Fax: 920-457-5773

### Contact Information:

Trang Adams / Regulatory Affairs Specialist  
1536 N. 18<sup>th</sup> Street  
Sheboygan, WI 53081  
Phone: 920-457-5051 Ext. 4251  
Fax: 920-457-5773  
E-Mail: tadams@americanortho.com

### Device Information:

Trade Name: Empower Clear  
Common Name: Ceramic Brackets  
Classification Name: Bracket, Ceramic, Orthodontic  
Classification Code: NJM  
Regulation Number: 872.5470

### Equivalent Legally Marketed Devices Information:

<u>510(k) #</u>	<u>Product Name</u>	<u>Device Manufacturer</u>
K080749	Radiance	American Orthodontics
K060837	In-Ovation C	Densply

### Description of the Device:

The Empower Clear line of products is single-use devices intended for use in conjunction with comprehensive orthodontics to control the movement of individual teeth. The Ceramic Bracket combines the aesthetics of a ceramic bracket with the versatility and ease of self ligation.

These brackets are comprised of several geometries that vary from bracket to bracket, corresponding to the intended tooth. These geometries contribute to the fit of the bracket to the tooth and also impart the axial control of the energy from the archwire.

### Indications for Use:

Empower Clear brackets are intended for orthodontic movement of teeth as diagnosed by an orthodontist. It is used temporarily and is removed upon completion of orthodontic treatment. Empower Clear brackets are intended to be single use only.



**Technological Characteristics Information:**

The material, Alumina Oxide [Al<sub>2</sub>O<sub>3</sub>], is used in the manufacturing of the Radiance brackets and In-Ovation C brackets – which is the same material used for Empower Clear.

The In-Ovation C clip material is the same as the Empower Clear clip material. The Radiance bracket does not have a clip.

The function and performance of the Empower Clear brackets are substantially equivalent to the predicate devices, as outlined in the following table.

Product Parameter	Device Name / Manufacturer			Substantial Equivalence Analysis
	Radiance / American Orthodontics	In-Ovation / Dentsply	Empower / American Orthodontics	
510(k) Number	K080749	K060837	Pending	N/A
Material	Al <sub>2</sub> O <sub>3</sub>	Al <sub>2</sub> O <sub>3</sub>	Al <sub>2</sub> O <sub>3</sub>	Equivalent
Intended Use	Orthodontic treatment is used to correct dental deficiencies and to improve the appearance of the patient. The brackets, arch wire and elastic o-rings form a force system that is designed to gradually move teeth into a normal alignment.	The Innovation C is intended for orthodontic movement of natural teeth, excluding the mandibular bicuspid teeth.	Ceramic Brackets are intended for orthodontic movement of teeth as diagnosed by an orthodontist. It is used temporarily and is removed upon completion of orthodontic treatment. Ceramic Brackets are intended to be single use only.	Equivalent
Single Use	YES	YES	YES	Equivalent
Non-Sterile Packaging	YES	YES	YES	Equivalent

**Biocompatibility Testing:**

Biocompatibility testing conducted on Empower Clear brackets indicates that the bracket material is safe for use. The bracket material, Aluminum Oxide [Al<sub>2</sub>O<sub>3</sub>], was found to be free of harmful extractables. No oral mucosa irritation or skin sensitization was detected with the material.

**Summary:**

The function and performance of Empower Clear bracket is similar to the predicates. There are no changes in the intended use and fundamental scientific technology. All of the materials used in the device have been used in legally marketed American Orthodontics devices. Minor differences in



technological characteristics do not raise new types of safety and effectiveness questions.



## DEPARTMENT OF HEALTH &amp; HUMAN SERVICES

## Public Health Service

U.S. Food and Drug Administration  
Center for Devices and Radiological Health  
Document Control Center WO66-G609  
10903 New Hampshire Avenue  
Silver Spring, MD 20993-0002

October 26, 2012

AMERICAN ORTHODONTICS  
1536 N. 18TH ST.  
SHEBOYGAN, WISCONSIN 53081  
ATTN: TRANG ADAMS

510k Number: K122753

Product: EMPOWER CLEAR

The additional information you have submitted has been received.

We will notify you when the processing of this submission has been completed or if any additional information is required. Please remember that all correspondence concerning your submission **MUST** be sent to the Document Mail Center at the above letterhead address. Correspondence sent to any address other than the one above will not be considered as part of your official premarket notification submission. Also, please note the new Blue Book Memorandum regarding Fax and E-mail Policy entitled, "Fax and E-Mail Communication with Industry about Premarket Files Under Review. Please refer to this guidance for information on current fax and e-mail practices at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm089402.htm>. On August 12, 2005 CDRH issued the Guidance for Industry and FDA Staff: Format for Traditional and Abbreviated 510(k)s. This guidance can be found at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm084365.htm>. Please refer to this guidance for assistance on how to format an original submission for a Traditional or Abbreviated 510(k).

The Safe Medical Devices Act of 1990, signed on November 28, states that you may not place this device into commercial distribution until you receive a letter from FDA allowing you to do so. As in the past, we intend to complete our review as quickly as possible. Generally we do so in 90 days. However, the complexity of a submission or a requirement for additional information may occasionally cause the review to extend beyond 90 days. Thus, if you have not received a written decision or been contacted within 90 days of our receipt date you may want to check with FDA to determine the status of your submission.

**Please ensure that whether you submit a 510(k) Summary as per 21 CFR 807.92, or a 510(k) Statement as per 21 CFR 807.93, it meets the content and format regulatory requirements.**

If you have procedural questions, please contact the Division of Small Manufacturers International and Consumer Assistance (DSMICA) at (301)796-7100 or at their toll-free number (800)638-2041, or contact the 510k staff at (301)796-5640.

Sincerely,

510(k) Staff

**Pugh, Dominique \***

---

**From:** Microsoft Outlook  
**To:** tadams@americanortho.com  
**Sent:** Friday, October 26, 2012 11:10 AM  
**Subject:** Relayed: K122753 AI Letter

**Delivery to these recipients or groups is complete, but no delivery notification was sent by the destination server:**

tadams@americanortho.com (tadams@americanortho.com)

Subject: K122753 AI Letter



October 23, 2012

12-13

K122753/S001

FDA CDRH DMC


OCT 26 2012

Received

RE: Empower Clear 510(k) # K122753

Please find enclosed two copies of the additional information requested sent on 10/22/12 via e-mail.

Sincerely,



Trang Adams / Regulatory Affairs Specialist

American Orthodontics

Phone: 920-457-5051 Ext. 4251

E-Mail: [tadams@americanortho.com](mailto:tadams@americanortho.com)

**EMPOWER CLEAR: 510(K) NUMBER K122753**

**ADDENDUM TO THE SUBMITTAL OF SEPTEMBER 12, 2012**

**CHEMICAL COMPOSITION:**

(b) (4)



**SYSTEM DESCRIPTION:**

The Empower Bracket line consists of Empower Self-Ligating Stainless Steel Brackets, Empower Reconvertible Stainless Steel Buccal Tubes and Empower Clear Brackets.

The Empower Clear is different from the other two Empower bracket lines, as it is not made of Stainless Steel, but is made of aluminum oxide material. Empower Clear is available as interactive or passive self-ligating systems. The interactive self-ligating system offers less friction at the beginning of treatment and more control in the later stages. The passive self-ligating system offers the benefits of less friction and lower force mechanics throughout treatment; no wires are actively engage by the clip.

The Empower Clear bracket system consists of 20 brackets: a group of 4 identical brackets (upper bicuspid), a second group of 4 identical brackets (lower centrals and laterals) and 12 unique brackets. Slot sizes available are 0.018" and 0.022" slot. The cuspid and bicuspid brackets are available with and without hooks that are used for auxillary attachments. The attached brochure lists three axis force system [torque, rotation, angulation]; as well as describing the identification scheme used for the system.

**BRACKET REMOVAL INSTRUCTIONS:**

Please refer to the proposed Empower Clear bonding and debonding (removal) instructions enclosed.

**BIOCOMPATIBILITY TESTING:**

(b) (4)




**COMPARISON OF RADIANCE & EMPOWER CLEAR BRACKETS:**

Similarities:

- 
- 
- 
- 
- 


(b) (4)



Differences:

- 
- 
- 
- 

(b)(4)





Empower Clear combines the aesthetics of a ceramic bracket with the versatility and ease of self ligation. Strong injection molded ceramic and aesthetically plated clips blend well with teeth, giving your patients the beautiful smile they deserve.

## A McLaughlin, Bennett, Trevisi System

INTERACTIVE

Maxillary	Torq	Ang	R/L	.018			.022		
				No Hook	Single Hook	Double Hook	No Hook	Single Hook	Double Hook
Central	+17	+4	R	1475-18-1117				1485-22-1117	
			L	1475-18-2117			1485-22-2117		
Lateral	+10	+8	R	1475-18-1210				1485-22-1210	
			L	1475-18-2210			1485-22-2210		
Cuspid	0	+8	R		1475-18-130HD	1475-18-130HU		1485-22-130HD	1485-22-130HU
			L		1475-18-230HD	1475-18-230HU		1485-22-230HD	1485-22-230HU
	-7	+8	R		1475-18-1307HD	1475-18-1307HU		1485-22-1307HD	1485-22-1307HU
			L		1475-18-2307HD	1475-18-2307HU		1485-22-2307HD	1485-22-2307HU
1st and 2nd Bicuspid	-7	0	R	1475-18-1407	1475-18-1407HD	1475-18-1407HU	1485-22-1407	1485-22-1407HD	1485-22-1407HU
			L	1475-18-2407	1475-18-2407HD	1475-18-2407HU	1485-22-2407	1485-22-2407HD	1485-22-2407HU



## A Roth System

INTERACTIVE

Maxillary	Torq	Ang	R/L	.018			.022		
				No Hook	Single Hook	Double Hook	No Hook	Single Hook	Double Hook
Central	+12	+5	R	1475-18-1112				1485-22-1112	
			L	1475-18-2112			1485-22-2112		
Lateral	+8	+9	R	1475-18-128				1485-22-128	
			L	1475-18-228			1485-22-228		
Cuspid	-2	+9	R		1475-18-1302HD	1475-18-1302HU		1485-22-1302HD	1485-22-1302HU
			L		1475-18-2302HD	1475-18-2302HU		1485-22-2302HD	1485-22-2302HU
1st and 2nd Bicuspid	-7	0	R	1475-18-1407	1475-18-1407HD	1475-18-1407HU	1485-22-1407	1485-22-1407HD	1485-22-1407HU
			L	1475-18-2407	1475-18-2407HD	1475-18-2407HU	1485-22-2407	1485-22-2407HD	1485-22-2407HU

### Bracket identification

Individual color coding ensures easy tooth identification and bonding accuracy



Color coded elastomeric Visual Positioning Aids in arch wire slots assist with bracket alignment during bonding. VPA's contain no latex and are easily removed with a scaler.

\*The American Orthodontics version of the McLaughlin, Bennett, Trevisi System is not claimed to be a duplication of any other, nor does American Orthodontics imply that it is endorsed in any way by Drs. McLaughlin, Bennett, or Trevisi.

\*The American Orthodontics version of the Roth System is not claimed to be a duplication of any other, nor does American Orthodontics imply that it is endorsed in any way by Dr. Roth.

# Empower Clear Braces

Empower Clear's patented Quad Matte™ bonding base provides a strong mechanical bond in the center of the base for reliable bond strength while delivering a weaker bond around the perimeter of the base for predictable debonding.



## Empower Clear Bonding

Due to the mechanical bond, no ceramic-specific bonding technique or adhesive is needed or recommended. Follow the adhesive manufacturer's recommendations for bonding of mechanical lock bases such as metal brackets.

### Direct Bonding

- > Bond Empower Clear brackets with the clip open.
- > Align the mesial and distal edges of the bracket base parallel to the long axis of the tooth.
- > Use the elastomeric Visual Positioning Aid (VPA) in the slot to align the bracket parallel to the desired occlusal plane.
- > If using light curable adhesive, cure Empower Clear through the face of the bracket for the amount of time recommended by the curing light manufacturer.
- > Remove the VPA using a scaler or other suitable instrument.



### Indirect Bonding

- > Remove the VPA from the slot and position bracket on the model with the clip closed.
- > Align the mesial and distal edges of the bracket base parallel to the long axis of the tooth.
- > Use the closed clip as a visual reference to align the bracket parallel to the desired occlusal plane.
- > To prevent impression material from entering the clip, fill the facial hole with lip balm or petroleum jelly.

### Empower Clear Rebonding

Always use a new bracket if repositioning or if replacing a premature debond.

### Empower Clear Debonding

- > Leave arch wires in place during debonding.
- > Remove adhesive flash from the mesial and distal edges of the bonding base.
- > Place blades of recommended debonding instrument along the enamel/adhesive interface, occlusal/gingival for upper centrals, mesial/distal for all other brackets.
- > Squeeze handles slowly until bracket releases from tooth.



001-301E  
Ceramic Bracket  
Debonding Instrument



AMERICAN  
ORTHODONTICS

3524 Washington Avenue Sheboygan, WI USA 53081 920-457-5051 USA and Canada: 1-800-558-7687  
info@americanortho.com www.americanortho.com © 2012 American Orthodontics Corporation

**Proudly made in the USA**

UTSL-34

# Empower<sup>®</sup> Self Ligating

Easy open and close.

Double-ended opening instrument provides two opening methods and works with the entire Empower family - Clear, metal, and molars.



## Gingival Opening Method

Single prong tip engages clip tail




## Facial Opening Method

Double prong tip engages facial hole\*




## Empower Tubes

Double prong tip engages facial hole\*



## Closing Clip

Clips close easily with just finger pressure



\* Use care when using the facial hole because over-opening can deform the clip. Use a finger to support the tooth end to prevent the clip from over-opening occlusally.

**Jones, Ashlee \***

---

**From:** Microsoft Outlook  
**To:** 'tadams@americanortho.com'  
**Sent:** Monday, October 22, 2012 9:57 AM  
**Subject:** Relayed: K122753 Hold Letter

Delivery to these recipients or groups is complete, but no delivery notification was sent by the destination server:

'tadams@americanortho.com' (tadams@americanortho.com) <mailto:tadams@americanortho.com>

Subject: K122753 Hold Letter



## DEPARTMENT OF HEALTH &amp; HUMAN SERVICES

## Public Health Service

U.S. Food and Drug Administration  
Center for Devices and Radiological Health  
Document Control Center WO66-G609  
10903 New Hampshire Avenue  
Silver Spring, MD 20993-0002

October 22, 2012

AMERICAN ORTHODONTICS  
1536 N. 18TH ST.  
SHEBOYGAN, WISCONSIN 53081  
ATTN: TRANG ADAMS

510k Number: K122753  
Product: EMPOWER CLEAR  
On Hold As of 10/19/2012

We are holding your above-referenced Premarket Notification (510(k)) for 30 days pending receipt of the additional information that was requested by the Office of Device Evaluation. Please remember that all correspondence concerning your submission MUST cite your 510(k) number and be sent in duplicate to the Document Mail Center at the above letterhead address. Correspondence sent to any address other than the one above will not be considered as part of your official premarket notification submission. Also, please note the new Blue Book Memorandum regarding Fax and E-mail Policy entitled, "Fax and E-Mail Communication with Industry about Premarket Files Under Review. Please refer to this guidance for information on current fax and e-mail practices at

<http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm089402.htm>.

The deficiencies identified represent the issues that we believe need to be resolved before our review of your 510(k) submission can be successfully completed. In developing the deficiencies, we carefully considered the statutory criteria as defined in Section 513(i) of the Federal Food, Drug, and Cosmetic Act for determining substantial equivalence of your device. We also considered the burden that may be incurred in your attempt to respond to the deficiencies. We believe that we have considered the least burdensome approach to resolving these issues. If, however, you believe that information is being requested that is not relevant to the regulatory decision or that there is a less burdensome way to resolve the issues, you should follow the procedures outlined in the "A Suggested Approach to Resolving Least Burdensome Issues" document. It is available on our Center web page at: <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/Overview/MedicalDeviceProvisionsofFDAModer nizationAct/ucm136685.htm>.

If after 30 days the additional information (AI), or a request for an extension of time, is not received, we will discontinue review of your submission and proceed to delete your file from our review system (21 CFR 807.87(l)). Please note our guidance document entitled, "Guidance for Industry and FDA Staff, FDA and Industry Actions on Premarket Notification (510(k)) Submissions: Effect on FDA Review Clock and Performance Assessment". If the submitter does submit a written request for an extension, FDA will permit the 510(k) to remain on hold for up to a maximum of 180 days from the date of the AI request. The purpose of this document is to assist agency staff and the device industry in understanding how various FDA and industry actions that may be taken on 510(k)s should affect the review clock for purposes of meeting the Medical Device User Fee and Modernization Act. You may review this document at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm089735.htm>. Pursuant to 21 CFR 20.29, a copy of your 510(k) submission will remain in the Office of Device Evaluation. If you then wish to resubmit this 510(k) notification, a new number will be assigned and your submission will be considered a new premarket notification submission.



Please remember that the Safe Medical Devices Act of 1990 states that you may not place this device into commercial distribution until you receive a decision letter from FDA allowing you to do so.

If you have procedural questions, please contact the Division of Small Manufacturers International and Consumer Assistance (DSMICA) at (301)796-7100 or at their toll-free number (800)638-2041, or contact the 510k staff at (301)796-5640.

Sincerely yours,

Marjorie Shulman  
Director, 510(k) Program  
Premarket Notification Section  
Office of Device Evaluation  
Center for Devices and Radiological Health



**COVER SHEET MEMORANDUM**

**From:** Reviewer Name Mary E Bue  
**Subject:** 510(k) Number K12A753  
**To:** The Record

Please list CTS decision code JH

- Refused to accept (Note: this is considered the first review cycle, See Screening Checklist [http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0\\_5631/Screening%20Checklist%207%202%2007.doc](http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0_5631/Screening%20Checklist%207%202%2007.doc))
- Hold (Additional Information or Telephone Hold)
- Final Decision (SE, SE with Limitations, NSE (select code below), Withdrawn, etc.).

**Not Substantially Equivalent (NSE) Codes**

- NO NSE for lack of predicate
- NI NSE for new intended use
- NQ NSE for new technology that raises new questions of safety and effectiveness
- NU NSE for new intended use AND new technology raising new questions of safety and effectiveness
- NP NSE for lack of performance data
- NS NSE no response
- NL NSE for lack of performance data AND no response
- NM NSE pre-amendment device call for PMAs (515i)
- NC NSE post-amendment device requires PMAs
- NH NSE for new molecular entity requires PMA
- TR NSE for transitional device

Please complete the following for a final clearance decision (i.e., SE, SE with Limitations, etc.):		YES	NO
Indications for Use Page	Attach IFU		
510(k) Summary /510(k) Statement	Attach Summary		
Truthful and Accurate Statement.	Must be present for a Final Decision		
Is the device Class III?			
If yes, does firm include Class III Summary?	Must be present for a Final Decision		
Does firm reference standards? (If yes, please attach form from <a href="http://www.fda.gov/opacom/morechoices/fdaforms/FDA-3654.pdf">http://www.fda.gov/opacom/morechoices/fdaforms/FDA-3654.pdf</a> )			
Is this a combination product? (Please specify category _____, see <a href="http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0_413b/COMBINATION%20PRODUCT%20ALGORITHM%20(REVISED%203-12-03).DOC">http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0_413b/COMBINATION%20PRODUCT%20ALGORITHM%20(REVISED%203-12-03).DOC</a> )			
Is this a reprocessed single use device? (Guidance for Industry and FDA Staff – MDUFMA - Validation Data in 510(k)s for Reprocessed Single-Use Medical Devices, <a href="http://www.fda.gov/cdrh/ode/guidance/1216.html">http://www.fda.gov/cdrh/ode/guidance/1216.html</a> )			
Is this device intended for pediatric use only?			
Is this a prescription device? (If both prescription & OTC, check both boxes.)			
Did the application include a completed FORM FDA 3674, Certification with Requirements of ClinicalTrials.gov Data Bank?			
Is clinical data necessary to support the review of this 510(k)?			
For United States-based clinical studies only: Did the application include a completed FORM FDA 3674, Certification with Requirements of ClinicalTrials.gov Data Bank? (If study was			

conducted in the United States, and FORM FDA 3674 was not included or incomplete, then applicant must be contacted to obtain completed form.)

Does this device include an Animal Tissue Source?

All Pediatric Patients age <= 21

Neonate/Newborn (Birth to 28 days)

Infant (29 days - < 2 years old)

Child (2 years - < 12 years old)

Adolescent (12 years - < 18 years old)

Transitional Adolescent A (18 - < 21 years old) Special considerations are being given to this group, different from adults age ≥ 21 (different device design or testing, different protocol procedures, etc.)

Transitional Adolescent B (18 - <= 21; No special considerations compared to adults => 21 years old)

Nanotechnology

Is this device subject to the Tracking Regulation? (Medical Device Tracking Guidance, <http://www.fda.gov/cdrh/comp/guidance/169.html>) Contact OC.

Regulation Number	Class*	Product Code
-------------------	--------	--------------

(\*If unclassified, see 510(k) Staff)

Additional Product Codes: \_\_\_\_\_

Review: \_\_\_\_\_

*Swan Russo*  
(Branch Chief)

*DSDB*  
(Branch Code)


*10/19/12*  
(Date)

Final Review: \_\_\_\_\_

*Swan Russo*  
(Division Director)

*10/19/12*  
(Date)


### 510(k) "SUBSTANTIAL EQUIVALENCE" DECISION-MAKING PROCESS



\* 510(k) Submissions compare new devices to marketed devices. FDA requests additional information if the relationship between marketed and "predicate" (pre-Amendments or reclassified post-Amendments) devices is unclear.

\*\* This decision is normally based on descriptive information alone, but limited testing information is sometimes required.

\*\*\* Data may be in the 510(k), other 510(k)s, the Center's classification files, or the literature.



To: The Record

From: Biologist, DEDB, DAGID, ODE, CDRH

Subject: K122753

Date: October 19, 2012

Background


American Orthodontics submitted a 510(k) for Empower Clear ceramic orthodontic brackets. I spoke with Ms. Trang Adams of AO today to request the following additional information:


1. the complete chemical composition of Empower Clear brackets;
2. describe the differences between the different Empower Clear brackets in the Empower line;
3. explain the recommended removal procedures for the brackets;
4. discuss the biocompatibility testing conducted for these brackets and include a summary of the testing in the 510(k) summary; and
5. discuss the differences between Empower Clear and the appropriate predicate device, and include a table with the physical properties of each device.

I explained to Ms. Adams that I will place her device on telephone hold until she submits this requested information to the Document Mail Center.

Recommendation

Place document on telephone hold.

  
Myra E. Browne

  
Browne 10/14/12



## DEPARTMENT OF HEALTH &amp; HUMAN SERVICES

## Public Health Service

U.S. Food and Drug Administration  
 Center for Devices and Radiological Health  
 Document Control Center WO66-G609  
 10903 New Hampshire Avenue  
 Silver Spring, MD 20993-0002

September 14, 2012

AMERICAN ORTHODONTICS  
 1536 N. 18TH ST.  
 SHEBOYGAN, WISCONSIN 53081  
 ATTN: TRANG ADAMS

510k Number: K122753

Received: 9/12/2012

Product: EMPOWER CLEAR

The Food and Drug Administration (FDA), Center for Devices and Radiological Health (CDRH), has received the Premarket Notification, (510(k)), you submitted in accordance with Section 510(k) of the Federal Food, Drug, and Cosmetic Act (Act) for the above referenced product and for the above referenced 510(k) submitter. Please note, if the 510(k) submitter is incorrect, please notify the 510(k) Staff immediately. We have assigned your submission a unique 510(k) number that is cited above. Please refer prominently to this 510(k) number in all future correspondence that relates to this submission. We will notify you when the processing of your 510(k) has been completed or if any additional information is required. **YOU MAY NOT PLACE THIS DEVICE INTO COMMERCIAL DISTRIBUTION UNTIL YOU RECEIVE A LETTER FROM FDA ALLOWING YOU TO DO SO.**

Please remember that all correspondence concerning your submission **MUST** be sent to the Document Mail Center (DMC) at the above letterhead address. Correspondence sent to any address other than the one above will not be considered as part of your official 510(k) submission.

On September 27, 2007, the President signed an act reauthorizing medical device user fees for fiscal years 2008 - 2012. The legislation - the Medical Device User Fee Amendments of 2007 is part of a larger bill, the Food and Drug Amendments Act of 2007. Please visit our website at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/Overview/MedicalDeviceUserFeeandModernizationActMDUFMA/default.htm>

for more information regarding fees and FDA review goals. In addition, effective January 2, 2008, any firm that chooses to use a standard in the review of ANY new 510(k) needs to fill out the new standards form (Form 3654) and submit it with their 510(k). The form may be found at <http://www.fda.gov/AboutFDA/ReportsManualsForms/Forms/default.htm>.

We remind you that Title VIII of the Food and Drug Administration Amendments Act of 2007 (FDAAA) amended the PHS Act by adding new section 402(j) (42 U.S.C. § 282(j)), which expanded the current database known as ClinicalTrials.gov to include mandatory registration and reporting of results for applicable clinical trials of human drugs (including biological products) and devices. Section 402(j) requires that a certification form <http://www.fda.gov/AboutFDA/ReportsManualsForms/Forms/default.htm> accompany 510(k)/HDE/PMA submissions. The agency has issued a draft guidance titled: "Certifications To Accompany Drug, Biological

Product, and Device Applications/Submissions: Compliance with Section 402(j) of The Public Health Service Act, Added By Title VIII of The Food and Drug Administration Amendments Act of 2007”

<http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/HowtoMarketYourDevice/PremarketSubmissions/PremarketNotification510k/ucm134034.htm>. According to the draft guidance, 510(k) submissions that do not contain clinical data do not need the certification form.

Please note the following documents as they relate to 510(k) review: 1) Guidance for Industry and FDA Staff entitled, “Interactive Review for Medical Device Submissions: 510(k)s, Original PMAs, PMA Supplements, Original BLAs and BLA Supplements”. This guidance can be found at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm089402.htm>. Please refer to this guidance for information on a formalized interactive review process. 2) Guidance for Industry and FDA Staff entitled, "Format for Traditional and Abbreviated 510(k)s". This guidance can be found at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm084365.htm>. Please refer to this guidance for assistance on how to format an original submission for a Traditional or Abbreviated 510(k).

In all future premarket submissions, we encourage you to provide an electronic copy of your submission. By doing so, you will save FDA resources and may help reviewers navigate through longer documents more easily. Under CDRH's e-Copy Program, you may replace one paper copy of any premarket submission (e.g., 510(k), IDE, PMA, HDE) with an electronic copy. For more information about the program, including the formatting requirements, please visit our web site at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/HowtoMarketYourDevice/PremarketSubmissions/ucm134508.html>. In addition, the 510(k) Program Video is now available for viewing on line at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/HowtoMarketYourDevice/PremarketSubmissions/PremarketNotification510k/ucm070201.htm>.

**Please ensure that whether you submit a 510(k) Summary as per 21 CFR 807.92, or a 510(k) Statement as per 21 CFR 807.93, it meets the content and format regulatory requirements.**

Lastly, you should be familiar with the regulatory requirements for medical devices available at Device Advice <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/default.htm>. If you have questions on the status of your submission, please contact DSMICA at (301)796-7100 or the toll-free number (800)638-2041, or at their internet address <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/default.htm>. If you have procedural questions, please contact the 510(k) Staff at (301)796-5640.

Sincerely,

510(k) Staff



**Williams, Michael \***

---

**From:** Microsoft Outlook  
**To:** tadams@americanortho.com  
**Sent:** Friday, September 14, 2012 8:40 AM  
**Subject:** Relayed: Ack Letter for K122753

**Delivery to these recipients or groups is complete, but no delivery notification was sent by the destination server:**

tadams@americanortho.com (tadams@americanortho.com)

Subject: Ack Letter for K122753



## DEPARTMENT OF HEALTH &amp; HUMAN SERVICES

Public Health Service

U.S. Food and Drug Administration  
 Center for Devices and Radiological Health  
 Document Control Center WO66-G609  
 10903 New Hampshire Avenue  
 Silver Spring, MD 20993-0002

September 10, 2012

USER FEE HOLD LETTER - HAVE NOT RECEIVED PAYM

AMERICAN ORTHODONTICS  
 1536 N. 18TH ST.  
 SHEBOYGAN, WISCONSIN 53081  
 ATTN: TRANG ADAMS

510k Number: K122753  
 Received: 9/7/2012  
 User Fee ID Number: 6063568  
 Product: EMPOWER CLEAR

The Food and Drug Administration (FDA) Center for Devices and Radiological Health (CDRH), has received the Premarket Notification you submitted in accordance with Section 510(k) of the Federal Food, Drug, and Cosmetic Act (Act) for the above referenced product. We have assigned your submission a unique 510(k) number that is cited above. Please refer prominently to this 510(k) number in all future correspondence that relates to this submission. **YOU MAY NOT PLACE THIS DEVICE INTO COMMERCIAL DISTRIBUTION UNTIL YOU RECEIVE A LETTER FROM FDA ALLOWING YOU TO DO SO.**

The Federal Food, Drug, and Cosmetic Act (the Act), as amended by the Medical Device User Fee and Modernization Act of 2002 (MDUFMA) and the FDA Amendments Act of 2007 (FDAAA) (Public Law 110-85), authorizes FDA to collect user fees for certain types of 510(k) submissions. The submission cannot be accepted for review until the fee is paid in full ; therefore, the file has been placed on hold. When your user fee payment has been received , review of the 510(k) will resume as of that date. Alternatively, you may request withdrawal of your submission. You now have the option to pay online by credit card. We recommend this form of payment. Credit card payments are directly linked to your user fee cover sheet and are processed the next business day. You may also pay by check. If you choose to mail a check, please send a check to one of the addresses listed below:

By Regular Mail  
 Food and Drug Administration  
 P.O. Box 956733  
 St. Louis, MO 63195-6733.

By Private Courier(e.g., Fed Ex, UPS, etc.)  
 U.S. Bank  
 956733  
 1005 Convention Plaza  
 St. Louis, MO 63101

The check should be made out to the Food and Drug Administration referencing the payment identification number, and a copy of the User Fee Cover sheet should be included with the check. A copy of the Medical Device User Fee Cover Sheet should be faxed to CDRH at (301)847-8120 referencing the 510(k) number if you have not already sent it in with your 510(k) submission. After the FDA has been notified of the receipt of your user fee payment, your 510(k) will be filed and the review will begin. If payment has not been received within 30 days, your 510(k) will be deleted from the system. Additional information on user fees and how to submit your user fee payment may be found at [www.fda.gov/cdrh/mdufma/fy09userfee.html](http://www.fda.gov/cdrh/mdufma/fy09userfee.html). In addition, the 510k Program Video is now available for viewing on line at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/HowtoMarketYourDevice/PremarketSubmissions/PremarketNotification510k/ucm070201.htm>.

In all future premarket submissions, we encourage you to provide an electronic copy of your submission. By doing so, you will save FDA resources and may help reviewers navigate through longer documents more easily. Under CDRH's e-Copy Program, you may replace one paper copy of any premarket submission (e.g., 510(k), IDE, PMA, or HDE) with an electronic copy. For more information about the program, including the formatting requirements, please visit our web site at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/HowtoMarketYourDevice/PremarketSubmissions/ucm134508.htm>.

Please note that since your 510(k) has not been reviewed, additional information may be required during the review process and the file may be placed on hold once again. If you are unsure as to whether or not you need to file a 510k Submission with FDA or what type of submission to submit, you should first telephone the Division of Small Manufacturers, International and Consumer Assistance (DSMICA), for guidance at (301)796-7100 or its toll-free number (800)638-2041, or contact them at their Internet address <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/default.htm>, or you may submit a 513(g) request for information regarding classification to the Document Mail Center at the address above. If you have any questions concerning receipt of your payment, please contact Edwena Jones at [Edwena.Jones@fda.hhs.gov](mailto:Edwena.Jones@fda.hhs.gov) or directly at (301)796-7200. If you have questions regarding the status of your 510(k) Submission, please contact DSMICA at the numbers or address above.

Sincerely yours,

Edwena Jones  
Consumer Safety Technician  
Premarket Notification Section  
Office of Device Evaluation  
Center for Devices and Radiological Health

**Mcdonald, Lisa \***

---

**From:** Microsoft Outlook  
**To:** tadams@americanortho.com  
**Sent:** Monday, September 10, 2012 10:24 AM  
**Subject:** Relayed: K122753 ACK Letter

**Delivery to these recipients or groups is complete, but no delivery notification was sent by the destination server:**

tadams@americanortho.com (tadams@americanortho.com)

Subject: K122753 ACK Letter

# American Orthodontics 510(k) Submission | 2012

---

## Table of Contents

Section 1: Medical Device User Fee .....	1
Medical Device User Fee Cover Sheet [MD# 6063568-956733] .....	1-1
Cover Sheet Order Confirmation [MD# 6063568-956733] .....	1-2
Section 2: CDRH Premarket Review Submission Cover Sheet.....	2
CDRH Premarket Review Submission Cover Sheet .....	2-1
Section 3: 510(k) Cover Letter .....	3
510(k) Cover Letter .....	3-1
Section 4: Indications for Use.....	4
Indications for Use Statement .....	4-1
Section 5: 510(k) Summary .....	5
510(k) Summary .....	5-1
Section 6: Truthful & Accuracy Statement .....	6
Truthful & Accuracy Statement.....	6-1
Section 7: Class III Summary & Certification .....	7
Not Applicable Statement .....	7-1
Section 8: Financial Certification or Disclosure Statement .....	8
Not Applicable Statement .....	8-1
Section 9: Declarations of Conformity & Summary Reports .....	9
Not Applicable Statement .....	9-1
Section 10: Executive Summary .....	10
Executive Summary .....	10-1
Section 11: Device Description.....	11
Device Description .....	11-1
Section 12: Substantial Equivalence Discussion.....	12
Substantial Equivalence Discussion .....	12-1

# American Orthodontics 510(k) Submission | 2012

---

## Table of Contents, Continued

Section 13: Proposed Labeling .....	13
Empower Clear Labels.....	13-1
Empower Clear Literature #1 .....	13-3
Empower Clear Literature #2 .....	13-4
Empower Clear Literature #3 .....	13-11
Section 14: Sterilization/ Shelf Life.....	14
Not Applicable Statement .....	14-1
Section 15: Biocompatibility.....	15
Biocompatibility .....	15-1
Section 16: Software .....	16
Not Applicable Statement .....	16-1
Section 17: Electromagnetic Compatibility / Electrical Safety .....	17
Not Applicable Statement .....	17-1
Section 18: Performance Testing - Bench .....	18
Not Applicable Statement .....	18-1
Section 19: Performance Testing - Animal .....	19
Not Applicable Statement .....	19-1
Section 20: Performance Testing - Clinical .....	20
Not Applicable Statement .....	20-1

**American Orthodontics 510(k) Submission | 2012**

---

**SECTION 1:  
MEDICAL DEVICE USER FEE**

American Orthodontics 510(k) Submission | 2012

---

2

**SECTION 2:**  
CDRH PREMARKET REVIEW  
SUBMISSION COVER SHEET



page 1-1

Form Approved: OMB No. 0910-511. See Instructions for OMB Statement.

DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION <b>MEDICAL DEVICE USER FEE COVER SHEET</b>		PAYMENT IDENTIFICATION NUMBER: (b)(4) Write the Payment Identification number on your check.
A completed cover sheet must accompany each original application or supplement subject to fees. If payment is sent by U.S. mail or courier, please include a copy of this completed form with payment. Payment and mailing instructions can be found at: <a href="http://www.fda.gov/oc/mdufma/coversheet.html">http://www.fda.gov/oc/mdufma/coversheet.html</a>		
1. COMPANY NAME AND ADDRESS (include name, street address, city state, country, and post office code)  AMERICAN ORTHODONTICS 1714 CAMBRIDGE AVE SHEBOYGAN WI 53082 US  1.1 EMPLOYER IDENTIFICATION NUMBER (EIN) *****3323	2. CONTACT NAME Trang Adams  2.1 E-MAIL ADDRESS tadams@americanortho.com  2.2 TELEPHONE NUMBER (include Area code) 920-457-5051 4251  2.3 FACSIMILE (FAX) NUMBER (Include Area code)	
3. TYPE OF PREMARKET APPLICATION (Select one of the following in each column; if you are unsure, please refer to the application descriptions at the following web site: <a href="http://www.fda.gov/oc/mdufma">http://www.fda.gov/oc/mdufma</a> ) <u>Select an application type:</u> <input checked="" type="checkbox"/> Premarket notification(510(k)); except for third party <input type="checkbox"/> 513(g) Request for Information <input type="checkbox"/> Biologics License Application (BLA) <input type="checkbox"/> Premarket Approval Application (PMA) <input type="checkbox"/> Modular PMA <input type="checkbox"/> Product Development Protocol (PDP) <input type="checkbox"/> Premarket Report (PMR) <input type="checkbox"/> Annual Fee for Periodic Reporting (APR) <input type="checkbox"/> 30-Day Notice		
3.1 Select a center <input checked="" type="checkbox"/> CDRH <input type="checkbox"/> CBER 3.2 Select one of the types below <input checked="" type="checkbox"/> Original Application <u>Supplement Types:</u> <input type="checkbox"/> Efficacy (BLA) <input type="checkbox"/> Panel Track (PMA, PMR, PDP) <input type="checkbox"/> Real-Time (PMA, PMR, PDP) <input type="checkbox"/> 180-day (PMA, PMR, PDP)		
4. ARE YOU A SMALL BUSINESS? (See the instructions for more information on determining this status) <input type="checkbox"/> YES, I meet the small business criteria and have submitted the required qualifying documents to FDA <input checked="" type="checkbox"/> NO, I am not a small business 4.1 If Yes, please enter your Small Business Decision Number:		
5. FDA WILL NOT ACCEPT YOUR SUBMISSION IF YOUR COMPANY HAS NOT PAID AN ESTABLISHMENT REGISTRATION FEE THAT IS DUE TO FDA. HAS YOUR COMPANY PAID ALL ESTABLISHMENT REGISTRATION FEES THAT ARE DUE TO FDA? <input checked="" type="checkbox"/> YES (All of our establishments have registered and paid the fee, or this is our first device, and we will register and pay the fee within 30 days of FDA's approval/clearance of this device.) <input type="checkbox"/> NO (If "NO," FDA will not accept your submission until you have paid all fees due to FDA. This submission will not be processed; see <a href="http://www.fda.gov/cdrh/mdufma">http://www.fda.gov/cdrh/mdufma</a> for additional information)		
6. IS THIS PREMARKET APPLICATION COVERED BY ANY OF THE FOLLOWING USER FEE EXCEPTIONS? IF SO, CHECK THE APPLICABLE EXCEPTION. <input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates <input type="checkbox"/> This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only <input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population <input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially		
7. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION FOR WHICH FEES WERE WAIVED DUE TO SOLE USE IN A PEDIATRIC POPULATION THAT NOW PROPOSES CONDITION OF USE FOR ANY ADULT POPULATION? (If so, the application is subject to the fee that applies for an original premarket approval application (PMA).) <input type="checkbox"/> YES <input checked="" type="checkbox"/> NO		
PAPERWORK REDUCTION ACT STATEMENT Public reporting burden for this collection of information is estimated to average 18 minutes per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to the address below.  Department of Health and Human Services, Food and Drug Administration, Office of Chief Information Officer, 1350 Piccard Drive, 4th Floor Rockville, MD 20850 [Please do NOT return this form to the above address, except as it pertains to comments on the burden estimate.]		
8. USER FEE PAYMENT AMOUNT SUBMITTED FOR THIS PREMARKET APPLICATION (b)(4)		27-Aug-2012

Form FDA 3601 (01/2007)

"Close Window" Print Cover sheet

Page 1-2



- [FAQ](#)
- [User Fees](#)
- [Draft Cover Sheet](#)
- [Previous Cover Sheets](#)
- [Profile](#)
- [Logout](#)

[Medical Device User Fee](#)

### Confirmation

YOUR PAYMENT IDENTIFICATION NUMBER IS: MD 6063568-956733

Your Cover Sheet has been submitted electronically. You must print and sign the hard copies. Include one in each copy of your application and include a copy with your payment.

Thank you for visiting the FDA User Fee Website. As part of our efforts to improve customer service, we would like to hear from you. Please [click here](#) to submit a survey. This will only take about 2 minutes to complete.

Coversheet	Creation Date	Last Update Date	
<a href="#">Medical Device User Fee and Modernization Act</a> (Print/View Final Coversheet)	1	27-AUG-2012	27-AUG-2012 Net: \$4,049.00

**Total: \$4,049.00**

#### Customer Information

Customer: AMERICAN ORTHODONTICS  
 Trang Adams  
 920-457-5051 4251  
 tadams@americanortho.com

#### Applicant Contact Information

Bill To: Trang Adams  
 AMERICAN ORTHODONTICS  
 1714 CAMBRIDGE AVE  
 SHEBOYGAN, WI 53082  
 UNITED STATES

[Pay Now](#) [Create Another Cover Sheet](#)

[Medical Device User Fee](#)  
[User Fees](#) | [Draft Cover Sheet](#) | [Previous Cover Sheets](#) | [Profile](#) | [Logout](#)

[FDA Home Page](#) | [Search FDA Site](#) | [FDA A-Z Index](#) | [Contact FDA](#) | [Privacy](#) | [Accessibility](#)

FDA Website Management Staff

*Trang Adams* 8/27/2012  
*S.W.M.* 8/27/12

Page 2-1

DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION <b>CDRH PREMARKET REVIEW SUBMISSION COVER SHEET</b>	Form Approval OMB No. 0910-0120 Expiration Date: December 31, 2013 See OMB Statement on page 5.
--	--

Date of Submission 09/04/2012	User Fee Payment ID Number (b)(4)	FDA Submission Document Number (if known) unknown
----------------------------------	--------------------------------------	--

**SECTION A TYPE OF SUBMISSION**

<b>PMA</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Premarket Report <input type="checkbox"/> Modular Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment <input type="checkbox"/> Licensing Agreement	<b>PMA &amp; HDE Supplement</b> <input type="checkbox"/> Regular (180 day) <input type="checkbox"/> Special <input type="checkbox"/> Panel Track (PMA Only) <input type="checkbox"/> 30-day Supplement <input type="checkbox"/> 30-day Notice <input type="checkbox"/> 135-day Supplement <input type="checkbox"/> Real-time Review <input type="checkbox"/> Amendment to PMA & HDE Supplement <input type="checkbox"/> Other	<b>PDP</b> <input type="checkbox"/> Original PDP <input type="checkbox"/> Notice of Completion <input type="checkbox"/> Amendment to PDP	<b>510(k)</b> <input checked="" type="checkbox"/> Original Submission: <input checked="" type="checkbox"/> Traditional <input type="checkbox"/> Special <input type="checkbox"/> Abbreviated (Complete section I, Page 5) <input type="checkbox"/> Additional Information <input type="checkbox"/> Third Party	<b>Meeting</b> <input type="checkbox"/> Pre-510(K) Meeting <input type="checkbox"/> Pre-IDE Meeting <input type="checkbox"/> Pre-PMA Meeting <input type="checkbox"/> Pre-PDP Meeting <input type="checkbox"/> Day 100 Meeting <input type="checkbox"/> Agreement Meeting <input type="checkbox"/> Determination Meeting <input type="checkbox"/> Other (specify):
<b>IDE</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement	<b>Humanitarian Device Exemption (HDE)</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment	<b>Class II Exemption Petition</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information	<b>Evaluation of Automatic Class III Designation (De Novo)</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information	<b>Other Submission</b> <input type="checkbox"/> 513(g) <input type="checkbox"/> Other (describe submission):

Have you used or cited Standards in your submission?  Yes  No (If Yes, please complete Section I, Page 5)

**SECTION B SUBMITTER, APPLICANT OR SPONSOR**

Company / Institution Name American Orthodontics		Establishment Registration Number (if known) 2126683	
Division Name (if applicable)		Phone Number (including area code) 920-457-5051	
Street Address 1536 N. 18th Street		FAX Number (including area code) 920-457-5773	
City Sheboygan	State / Province Wisconsin	ZIP/Postal Code 53081	Country U.S.A
Contact Name Trang Adams			
Contact Title Regulatory Affairs Specialist		Contact E-mail Address tadams@americanortho.com	

**SECTION C APPLICATION CORRESPONDENT (e.g., consultant, if different from above)**

Company / Institution Name			
Division Name (if applicable)		Phone Number (including area code)	
Street Address		FAX Number (including area code)	
City	State / Province	ZIP Code	Country
Contact Name			
Contact Title		Contact E-mail Address	

Page 2-2

SECTION D1			REASON FOR APPLICATION - PMA, PDP, OR HDE		
<input type="checkbox"/> New Device <input type="checkbox"/> Withdrawal <input type="checkbox"/> Additional or Expanded Indications <input type="checkbox"/> Request for Extension <input type="checkbox"/> Post-approval Study Protocol <input type="checkbox"/> Request for Applicant Hold <input type="checkbox"/> Request for Removal of Applicant Hold <input type="checkbox"/> Request to Remove or Add Manufacturing Site	<input type="checkbox"/> Change in design, component, or specification: <input type="checkbox"/> Software / Hardware <input type="checkbox"/> Color Additive <input type="checkbox"/> Material <input type="checkbox"/> Specifications <input type="checkbox"/> Other (specify below)	<input type="checkbox"/> Location change: <input type="checkbox"/> Manufacturer <input type="checkbox"/> Sterilizer <input type="checkbox"/> Packager			
<input type="checkbox"/> Process change: <input type="checkbox"/> Manufacturing <input type="checkbox"/> Packaging <input type="checkbox"/> Sterilization <input type="checkbox"/> Other (specify below)	<input type="checkbox"/> Labeling change: <input type="checkbox"/> Indications <input type="checkbox"/> Instructions <input type="checkbox"/> Performance Characteristics <input type="checkbox"/> Shelf Life <input type="checkbox"/> Trade Name <input type="checkbox"/> Other (specify below)	<input type="checkbox"/> Report Submission: <input type="checkbox"/> Annual or Periodic <input type="checkbox"/> Post-approval Study <input type="checkbox"/> Adverse Reaction <input type="checkbox"/> Device Defect <input type="checkbox"/> Amendment			
<input type="checkbox"/> Response to FDA correspondence:		<input type="checkbox"/> Change in Ownership <input type="checkbox"/> Change in Correspondent <input type="checkbox"/> Change of Applicant Address			
<input type="checkbox"/> Other Reason (specify):					

SECTION D2			REASON FOR APPLICATION - IDE		
<input type="checkbox"/> New Device <input type="checkbox"/> New Indication <input type="checkbox"/> Addition of Institution <input type="checkbox"/> Expansion / Extension of Study <input type="checkbox"/> IRB Certification <input type="checkbox"/> Termination of Study <input type="checkbox"/> Withdrawal of Application <input type="checkbox"/> Unanticipated Adverse Effect <input type="checkbox"/> Notification of Emergency Use <input type="checkbox"/> Compassionate Use Request <input type="checkbox"/> Treatment IDE <input type="checkbox"/> Continued Access	<input type="checkbox"/> Change in: <input type="checkbox"/> Correspondent / Applicant <input type="checkbox"/> Design / Device <input type="checkbox"/> Informed Consent <input type="checkbox"/> Manufacturer <input type="checkbox"/> Manufacturing Process <input type="checkbox"/> Protocol - Feasibility <input type="checkbox"/> Protocol - Other <input type="checkbox"/> Sponsor	<input type="checkbox"/> Response to FDA Letter Concerning: <input type="checkbox"/> Conditional Approval <input type="checkbox"/> Deemed Approved <input type="checkbox"/> Deficient Final Report <input type="checkbox"/> Deficient Progress Report <input type="checkbox"/> Deficient Investigator Report <input type="checkbox"/> Disapproval <input type="checkbox"/> Request Extension of Time to Respond to FDA <input type="checkbox"/> Request Meeting <input type="checkbox"/> Request Hearing			
<input type="checkbox"/> Report submission: <input type="checkbox"/> Current Investigator <input type="checkbox"/> Annual Progress Report <input type="checkbox"/> Site Waiver Report <input type="checkbox"/> Final					
<input type="checkbox"/> Other Reason (specify):					

SECTION D3			REASON FOR SUBMISSION - 510(k)		
<input checked="" type="checkbox"/> New Device	<input type="checkbox"/> Additional or Expanded Indications	<input type="checkbox"/> Change in Technology			
<input type="checkbox"/> Other Reason (specify):					

Page 2-3

**SECTION E ADDITIONAL INFORMATION ON 510(K) SUBMISSIONS**

Product codes of devices to which substantial equivalence is claimed				Summary of, or statement concerning, safety and effectiveness information <input checked="" type="checkbox"/> 510 (k) summary attached <input type="checkbox"/> 510 (k) statement
1 NJM	2	3	4	
5	6	7	8	

Information on devices to which substantial equivalence is claimed (if known)		
510(k) Number	Trade or Proprietary or Model Name	Manufacturer
1 K080749	1 Radiance	1 American Orthodontics
2 K060837	2 In-Ovation C	2 Dctnsply
3	3	3
4	4	4
5	5	5
6	6	6

**SECTION F PRODUCT INFORMATION - APPLICATION TO ALL APPLICATIONS**

Common or usual name or classification name  
 Ceramic Bracket

Trade or Proprietary or Model Name for This Device	Model Number
1 Empower Clear	1
2	2
3	3
4	4
5	5

FDA document numbers of all prior related submissions (regardless of outcome)					
1	2	3	4	5	6
7	8	9	10	11	12

Data Included in Submission  
 Laboratory Testing       Animal Trials       Human Trials

**SECTION G PRODUCT CLASSIFICATION - APPLICATION TO ALL APPLICATIONS**

Product Code NJM	C.F.R. Section (if applicable) 872.5470	Device Class <input type="checkbox"/> Class I <input checked="" type="checkbox"/> Class II <input type="checkbox"/> Class III <input type="checkbox"/> Unclassified
Classification Panel		

Indications (from labeling)  
 Indications for Use:  
 Empower Clear brackets are intended for orthodontic movement of teeth as diagnosed by an orthodontist. It is used temporarily and is removed upon completion of orthodontic treatment. Empower Clear brackets are intended to be single use only.

Page 2-4

<p><b>Note:</b> Submission of the information entered in Section H does not affect the need to submit device establishment registration.</p>		<p>FDA Document Number (if known)</p> <p>unknown</p>	
<p><b>SECTION H MANUFACTURING / PACKAGING / STERILIZATION SITES RELATING TO A SUBMISSION</b></p>			
<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete		Facility Establishment Identifier (FEI) Number	
		<input checked="" type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Contract Manufacturer <input type="checkbox"/> Repackager / Relabeler	
Company / Institution Name American Orthodontics		Establishment Registration Number 2126683	
Division Name (if applicable)		Phone Number (including area code) 920-457-5051	
Street Address 1536 N. 18th Street		FAX Number (including area code) 920-457-5773	
City Sheboygan		State / Province Wisconsin	ZIP Code 53081
		Country U.S.A.	
Contact Name Trang Adams		Contact Title Regulatory Affairs Specialist	
		Contact E-mail Address tadams@americanortho.com	
<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete		Facility Establishment Identifier (FEI) Number	
		<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Contract Manufacturer <input type="checkbox"/> Repackager / Relabeler	
Company / Institution Name		Establishment Registration Number	
Division Name (if applicable)		Phone Number (including area code)	
Street Address		FAX Number (including area code)	
City		State / Province	ZIP Code
		Country	
Contact Name		Contact Title	
		Contact E-mail Address	
<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete		Facility Establishment Identifier (FEI) Number	
		<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Contract Manufacturer <input type="checkbox"/> Repackager / Relabeler	
Company / Institution Name		Establishment Registration Number	
Division Name (if applicable)		Phone Number (including area code)	
Street Address		FAX Number (including area code)	
City		State / Province	ZIP Code
		Country	
Contact Name		Contact Title	
		Contact E-mail Address	

Page 2-5

**SECTION I UTILIZATION OF STANDARDS**

**Note:** Complete this section if your application or submission cites standards or includes a "Declaration of Conformity to a Recognized Standard" statement.

	Standards No.	Standards Organization	Standards Title	Version	Date
1					
2					
3					
4					
5					
6					
7					

**Please include any additional standards to be cited on a separate page.**

**Public reporting burden for this collection of information** is estimated to average 0.5 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to:

Department of Health and Human Services  
 Food and Drug Administration  
 Office of Chief Information Officer  
 1350 Piccard Drive, Room 400  
 Rockville, MD 20850

*An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number.*

American Orthodontics 510(k) Submission | 2012

---

**SECTION 3:**  
**510(K) COVER LETTER**



K122753



AMERICAN  
ORTHODONTICS

### 510(k) Cover Letter

**Submission Date: September 4, 2012**

To Whom It May Concern,

This 510(k) Submission Cover Letter is submitted for the following company:

**Company:**

American Orthodontics

**Establishment Registration Number:**

2126683

FDA CDRH DMC

SEP - 7 2012

Received *Kelo*

This 510(k) Submission is **Submitted By:**

Trang Adams / Regulatory Affairs Specialist  
1536 N. 18<sup>th</sup> Street  
Sheboygan, WI 53081  
Phone: 920-457-5051 Ext. 4251  
Fax: 920-457-5773

The **Common Name of the Device** is:

Ceramic Bracket

The **Trade Name of the Device** is:

Empower Clear

The **Classification Name of the Device** is:

Bracket, Ceramic, Orthodontic [NJM]

The **510(k) Submission Reason** is:

New Device

The following are **Legally Marketed Devices [Predicate]** that is considered Substantially Equivalent [SE]:

<u>510(k) #</u>	<u>Product Name</u>
K080749	Radiance
K060837	In-Ovation C

The **Compliance with Any Special Controls** are as follows:

No applicable mandatory performance standards or special controls exist for this device.

Sincerely,

Trang Adams / Regulatory Affairs Specialist

American Orthodontics 510(k) Submission | 2012

---

**SECTION 4:**  
**INDICATIONS FOR USE**

R 122753



## Indications for Use Statement

**510(k) Number (if known):**

Unknown

**Device Name:**

Empower Clear

**Indications for Use:**

Empower Clear brackets are intended for orthodontic movement of teeth as diagnosed by an orthodontist. It is used temporarily and is removed upon completion of orthodontic treatment. Empower Clear brackets are intended to be single use only.

**Prescription Use And/Or Over-The-Counter Use:**

- Prescription use by orthodontist only
- Not available Over-The-Counter [OTC]

American Orthodontics 510(k) Submission | 2012

---

**SECTION 5:**  
510(K) SUMMARY

5



## 510(k) Summary

Page 5-1

**Preparation Date: August 28, 2012**

### **Company Information:**

American Orthodontics  
1714 Cambridge Avenue  
Sheboygan, WI 53081  
Phone: 920-457-5051  
Fax: 920-457-5773

### **Contact Information:**

Trang Adams / Regulatory Affairs Specialist  
1536 N. 18<sup>th</sup> Street  
Sheboygan, WI 53081  
Phone: 920-457-5051 Ext. 4251  
Fax: 920-457-5773  
E-Mail: tadams@americanortho.com

### **Device Information:**

Trade Name: Empower Clear  
Common Name: Ceramic Brackets  
Classification Name: Bracket, Ceramic, Orthodontic  
Classification Code: NJM  
Regulation Number: 872.5470

### **Equivalent Legally Marketed Devices Information:**

<u>510(k) #</u>	<u>Product Name</u>	<u>Device Manufacturer</u>
K080749	Radiance	American Orthodontics
K060837	In-Ovation C	Densply

### **Description of the Device:**

The Empower Clear line of products is single-use devices intended for use in conjunction with comprehensive orthodontics to control the movement of individual teeth. The Ceramic Bracket combines the aesthetics of a ceramic bracket with the versatility and ease of self ligation.

These brackets are comprised of several geometries that vary from bracket to bracket, corresponding to the intended tooth. These geometries contribute to the fit of the bracket to the tooth and also impart the axial control of the energy from the archwire.

### **Indications for Use:**

Empower Clear brackets are intended for orthodontic movement of teeth as diagnosed by an orthodontist. It is used temporarily and is removed upon completion of orthodontic treatment. Empower Clear brackets are intended to be single use only.



**Technological Characteristics Information:**

The material, Alumina Oxide [Al<sub>2</sub>O<sub>3</sub>], is used in the manufacturing of the Radiance brackets and In-Ovation C brackets – which is the same material used for Empower Clear.

The In-Ovation C clip material is the same as the Empower Clear clip material. The Radiance bracket does not have a clip.

The function and performance of the Empower Clear brackets are substantially equivalent to the predicate devices, as outlined in the following table.

Product Parameter	Device Name / Manufacturer			Substantial Equivalence Analysis
	Radiance / American Orthodontics	In-Ovation / Dentsply	Empower / American Orthodontics	
510(k) Number	K080749	K060837	Pending	N/A
Material	Al <sub>2</sub> O <sub>3</sub>	Al <sub>2</sub> O <sub>3</sub>	Al <sub>2</sub> O <sub>3</sub>	Equivalent
Intended Use	Orthodontic treatment is used to correct dental deficiencies and to improve the appearance of the patient. The brackets, arch wire and elastic o-rings form a force system that is designed to gradually move teeth into a normal alignment.	The Innovation C is intended for orthodontic movement of natural teeth, excluding the mandibular bicuspid teeth.	Ceramic Brackets are intended for orthodontic movement of teeth as diagnosed by an orthodontist. It is used temporarily and is removed upon completion of orthodontic treatment. Ceramic Brackets are intended to be single use only.	Equivalent
Single Use	YES	YES	YES	Equivalent
Non-Sterile Packaging	YES	YES	YES	Equivalent

**Summary:**

The function and performance of Empower Clear bracket is similar to the predicates. There are no changes in the intended use and fundamental scientific technology. All of the materials used in the device have been used in legally marketed American Orthodontics devices. Minor differences in technological characteristics do not raise new types of safety and effectiveness questions.

**American Orthodontics 510(k) Submission | 2012**

---

**SECTION 6:  
TRUTHFUL & ACCURACY STATEMENT**

6



## Truthful & Accuracy Statement

Page 6-1

510(k) Number: Unknown

### Premarket Notification Truthful & Accurate Statement [As required by 21 CFR 807.87(k)]

I certify that, in my capacity as Regulatory Affairs Specialist of American Orthodontics, I believe to the best of my knowledge, that all data and information submitted in the premarket notification are truthful and accurate and that no material fact has been omitted.

*Trang Adams 8/28/2012*

Trang Adams / Regulatory Affairs Specialist



**American Orthodontics 510(k) Submission | 2012**

---

**SECTION 7:  
CLASS III SUMMARY & CERTIFICATION**



**Page 7-1**

**Class III Summary & Certification**

**\*\*This section is not applicable. \*\***

American Orthodontics 510(k) Submission | 2012

---

8

**SECTION 8:**  
**FINANCIAL CERTIFICATION OR**  
**DISCLOSURE STATEMENT**



**Page 8-1**

**Financial Certification or Disclosure Statement**

**\*\*This section is not applicable. \*\***

American Orthodontics 510(k) Submission | 2012

---

**SECTION 9:**  
**DECLARATIONS OF CONFORMITY &  
SUMMARY REPORTS**



**Page 9-1**

**Declaration of Conformity & Summary Reports**

**\*\*This section is not applicable. \*\***

American Orthodontics 510(k) Submission | 2012

---

**SECTION 10:**  
**EXECUTIVE SUMMARY**



## Executive Summary

### Device Information:

Trade Name: Empower Clear  
Common Name: Ceramic Brackets  
Classification Name: Bracket, Ceramic, Orthodontic  
Classification Code: NJM  
Regulation Number: 872.5470

### Description of the Device:

The Empower Clear line of products is single-use devices intended for use in conjunction with comprehensive orthodontics to control the movement of individual teeth. The Empower Clear combines the aesthetics of a ceramic bracket with the versatility and ease of self ligation.

These brackets are comprised of several geometries that vary from bracket to bracket, corresponding to the intended tooth. These geometries contribute to the fit of the bracket to the tooth and also impart the axial control of the energy from the archwire.

### Indications for Use:

Empower Clear brackets are intended for orthodontic movement of teeth as diagnosed by an orthodontist. It is used temporarily and is removed upon completion of orthodontic treatment. Empower Clear brackets are intended to be single use only.

### Summary:

The function and performance of Empower Clear bracket is similar to the predicates, as outlined below in "Technological Characteristics Information". There are no changes in the intended use and fundamental scientific technology. All of the materials used in the device have been used in legally marketed American Orthodontics devices. Minor differences in technological characteristics do not raise new types of safety and effectiveness questions.





**Technological Characteristics Information:**

The material, Alumina Oxide [Al<sub>2</sub>O<sub>3</sub>], is used in the manufacturing of the Radiance brackets and In-Ovation C brackets – which is the same material used for Empower Clear.

The In-Ovation C clip material is the same as the Empower Clear clip material. The Radiance bracket does not have a clip.

The function and performance of the Empower Clear brackets are substantially equivalent to the predicate devices, as outlined in the following table.

Product Parameter	Device Name / Manufacturer			Substantial Equivalence Analysis
	Radiance / American Orthodontics	In-Ovation / Dentsply	Empower / American Orthodontics	
510(k) Number	K080749	K060837	Pending	N/A
Material	Al <sub>2</sub> O <sub>3</sub>	Al <sub>2</sub> O <sub>3</sub>	Al <sub>2</sub> O <sub>3</sub>	Equivalent
Intended Use	Orthodontic treatment is used to correct dental deficiencies and to improve the appearance of the patient. The brackets, arch wire and elastic o-rings form a force system that is designed to gradually move teeth into a normal alignment.	The Innovation C is intended for orthodontic movement of natural teeth, excluding the mandibular bicuspid teeth.	Ceramic Brackets are intended for orthodontic movement of teeth as diagnosed by an orthodontist. It is used temporarily and is removed upon completion of orthodontic treatment. Ceramic Brackets are intended to be single use only.	Equivalent
Single Use	YES	YES	YES	Equivalent
Non-Sterile Packaging	YES	YES	YES	Equivalent

American Orthodontics 510(k) Submission | 2012

---

**SECTION 11:**  
**DEVICE DESCRIPTION**



## Device Description

Page 11-1

### Device Information:

Trade Name: Empower Clear  
Common Name: Ceramic Brackets  
Classification Name: Bracket, Ceramic, Orthodontic  
Classification Code: NJM  
Regulation Number: 872.5470

### Description of the Device:

The Empower Clear bracket line of products is single-use devices intended for use in conjunction with comprehensive orthodontics to control the movement of individual teeth. The Ceramic Bracket combines the aesthetics of a ceramic bracket with the versatility and ease of self ligation.

These brackets are comprised of several geometries that vary from bracket to bracket, corresponding to the intended tooth. These geometries contribute to the fit of the bracket to the tooth and also impart the axial control of the energy from the archwire.

(b) (4)



(b)(4) Engineering Drawings



(b)(4) Engineering Drawings



(b)(4) Engineering Drawings



(b)(4) Engineering Drawings



(b)(4) Engineering Drawings





(b)(4) Engineering Drawings



(b)(4) Engineering Drawings



American Orthodontics 510(k) Submission | 2012

---

**SECTION 12:**  
SUBSTANTIAL EQUIVALENCE  
DISCUSSION

12



## Substantial Equivalence Discussion

Page 12-1

### Device Information:

Trade Name: Empower Clear  
 Common Name: Ceramic Brackets  
 Classification Name: Bracket, Ceramic, Orthodontic  
 Classification Code: NJM  
 Regulation Number: 872.5470

### Indications for Use:

Empower Clear brackets are intended for orthodontic movement of teeth as diagnosed by an orthodontist. It is used temporarily and is removed upon completion of orthodontic treatment. Empower Clear brackets are intended to be single use only.

### Equivalent Legally Marketed Devices Information:

<u>510(k) #</u>	<u>Product Name</u>	<u>Device Manufacturer</u>
K080749	Radiance	American Orthodontics
K060837	In-Ovation C	Densply

### Technological Characteristics Information:

The material, Alumina Oxide [Al<sub>2</sub>O<sub>3</sub>], is used in the manufacturing of the Radiance brackets and In-Ovation C brackets – which is the same material used for Empower Clear.

The In-Ovation C clip material is the same as the Empower Clear clip material. The Radiance bracket does not have a clip.

The function and performance of the Empower Clear brackets are substantially equivalent to the predicate devices, as outlined in the following table.

### Summary:

The function and performance of Empower Clear bracket is similar to the predicates. There are no changes in the intended use and fundamental scientific technology. All of the materials used in the device have been used in legally marketed American Orthodontics devices. Minor differences in technological characteristics do not raise new types of safety and effectiveness questions.



Technological characteristic comparisons are outlined below:

Product Parameter	Device Name / Manufacturer			Substantial Equivalence Analysis
	Radiance / American Orthodontics	In-Ovation / Dentsply	Empower / American Orthodontics	
510(k) Number	K080749	K060837	Pending	N/A
Material	Al <sub>2</sub> O <sub>3</sub>	Al <sub>2</sub> O <sub>3</sub>	Al <sub>2</sub> O <sub>3</sub>	Equivalent
Intended Use	Orthodontic treatment is used to correct dental deficiencies and to improve the appearance of the patient. The brackets, arch wire and elastic o-rings form a force system that is designed to gradually move teeth into a normal alignment.	The Innovation C is intended for orthodontic movement of natural teeth, excluding the mandibular bicuspid teeth.	Ceramic Brackets are intended for orthodontic movement of teeth as diagnosed by an orthodontist. It is used temporarily and is removed upon completion of orthodontic treatment. Ceramic Brackets are intended to be single use only.	Equivalent
Single Use	YES	YES	YES	Equivalent
Non-Sterile Packaging	YES	YES	YES	Equivalent

**American Orthodontics 510(k) Submission | 2012**

---

**SECTION 13:  
PROPOSED LABELING**

13

**AO** AMERICAN  
ORTHODONTICS

## Empower<sup>®</sup> Self Ligating System

.018     .022

Interactive

Passive

CE 0408 Contains Nickel and/or Chromium Patents #5,630,715 #5,722,826 #4,659,309

EC REP EC Certification Service GmbH • Sandgasse 39a, A-9300 St. Veit/Glan, Austria  
Made in USA

American Orthodontics • 1714 Cambridge Avenue • Sheboygan, WI 53081 • 1-800-558-7687

 American Orthodontics • 1714 Cambridge Avenue • Sheboygan, WI USA 53081 • 1-800-558-7687




**Empower**  
Clear Braces

Rx/Patient \_\_\_\_\_

- .018     .022
- Interactive
- Passive

CE 0408

 EC Certification Service GmbH • Sandgasse39a, A-9300 • St. Veit/Glan, Austria    Made in USA



Contains Nickel and/or Chromium

Patents #5,630,715  
#5,722,826  
#4,659,309

L-000424  
REV0





Empower Clear combines the aesthetics of a ceramic bracket with the versatility and ease of self ligation. Strong injection molded ceramic and aesthetically plated clips blend well with teeth, giving your patients the beautiful smile they deserve.

### A McLaughlin, Bennett, Trevisi System



Maxillary	Torq	Ang	R/L	.018			.022		
				No Hook	Single Hook	Double Hook	No Hook	Single Hook	Double Hook
Central	+17	+4	R	1475-18-1117			1485-22-1117		
			L	1475-18-2117			1485-22-2117		
Lateral	+10	+8	R	1475-18-1210			1485-22-1210		
			L	1475-18-2210			1485-22-2210		
Cuspid	0	+8	R		1475-18-130HD	1475-18-130HU		1485-22-130HD	1485-22-130HU
			L		1475-18-230HD	1475-18-230HU		1485-22-230HD	1485-22-230HU
	-7	+8	R		1475-18-1307HD	1475-18-1307HU		1485-22-1307HD	1485-22-1307HU
			L		1475-18-2307HD	1475-18-2307HU		1485-22-2307HD	1485-22-2307HU
1st and 2nd Bicuspid	-7	0	R	1475-18-1407	1475-18-1407HD	1475-18-1407HU	1485-22-1407	1485-22-1407HD	1485-22-1407HU
			L	1475-18-2407	1475-18-2407HD	1475-18-2407HU	1485-22-2407	1485-22-2407HD	1485-22-2407HU



### A Roth System



Maxillary	Torq	Ang	R/L	.018			.022		
				No Hook	Single Hook	Double Hook	No Hook	Single Hook	Double Hook
Central	+12	+5	R	1475-18-1112			1485-22-1112		
			L	1475-18-2112			1485-22-2112		
Lateral	+8	+9	R	1475-18-128			1485-22-128		
			L	1475-18-228			1485-22-228		
Cuspid	-2	+9	R		1475-18-1302HD	1475-18-1302HU		1485-22-1302HD	1485-22-1302HU
			L		1475-18-2302HD	1475-18-2302HU		1485-22-2302HD	1485-22-2302HU
1st and 2nd Bicuspid	-7	0	R	1475-18-1407	1475-18-1407HD	1475-18-1407HU	1485-22-1407	1485-22-1407HD	1485-22-1407HU
			L	1475-18-2407	1475-18-2407HD	1475-18-2407HU	1485-22-2407	1485-22-2407HD	1485-22-2407HU

#### Bracket Identification


Individual color coding ensures easy tooth identification and bonding accuracy



Color coded elastomeric Visual Positioning Aids in arch wire slots assist with bracket alignment during bonding. VPA's contain no latex and are easily removed with a scaler.

\*The American Orthodontics version of the McLaughlin, Bennett, Trevisi System is not claimed to be a duplication of any other, nor does American Orthodontics imply that it is endorsed in any way by Drs. McLaughlin, Bennett, or Trevisi.

\*The American Orthodontics version of the Roth System is not claimed to be a duplication of any other, nor does American Orthodontics imply that it is endorsed in any way by Dr. Roth.



**Empower**  
Self Ligating




**AO**  
AMERICAN  
ORTHODONTICS


Page 13-4



## It's all about choice.

Empower is the industry's most complete self ligating bracket system, bringing you choice like you've never had before. Empower's unique design options not only adapt to your treatment philosophies, they allow treatment to be tailored to your patients' needs. This unprecedented array of possibilities inspires you to treat patients your way. Empower puts control back in your hands.

-  Clear bracket combines performance with patient-pleasing aesthetics
-  Metal brackets offer fully interactive, fully passive, or combination Dual Activation™ options
-  Reconvertible molar tubes provide versatility without compromise



### Quality Features

Empower's complete line of brackets - clear, metal, and molar tubes - gives you the tools you need to give your patients the smile they deserve.



- > Low profile, patient-friendly hooks (top oval) >
  - Setting: hygienic, symmetrical hooks and wires
  - Use: a desirable white, stain-fighting polymer coating
- > Compound contoured Quad-Mate™ base >
  - Allows strong bands to connect and wear on edges for easy and predictable bonding
- > Chamfered slot entrances >
  - Eliminates sharp corners to reduce wire bending
- > Durable, aesthetically plated clip >
  - Black-velvet strands = 100% enamel for easier and smoother placement
- > Individual color coding and Visual Positioning Aids >
  - Facilitates easy tooth placement and bonding strategy
- < Clip tail >
  - Provides a secure original bonding strand
- < Smooth rounded contours >
  - Delivers exceptional patient comfort
- < Facial hole >
  - Adds local bonding method when clip tail is inaccessible
- < Interactive clip design >
  - Contains lower friction mechanics early in treatment with increased control during the critical finishing phase
- < Ceramic Injection Molded bracket body >
  - Results in a strong, very aesthetic bracket



- > Low profile, patient-friendly hooks (top oval) >
  - Delivers the hook design you need to hold and load every case with comfort
- > Compound contoured Maximum Retention™ pads >
  - Compares reliable bond strength with predictable bonding
- > Rounded slot entrances >
  - Eliminates sharp corners to reduce wire bending
- > Extended tie wings >
  - Provides easy flexibility and allows differential spring
- < Center Line >
  - Allows accurate bonding alignment with long axis of tooth
- < Keypast notch >
  - Facets wire engaged in a slot even with heavily rotated or maloccluded
- < Precision heat treated chromium cobalt clip >
  - Provides dependable performance throughout treatment
- < Individual color coding >
  - Facilitates easy tooth identification during treatment



- > Multiple design choices >
  - Provides single and double tube options, in no band or well style
- > Funneled tube entrances >
  - Reduces wire bending and allows easier wire insertion
- > Easy access facial hole >
  - Simplifies opening the clip significantly
- < Low profile, patient-friendly hooks >
  - Provides maximum comfort and ease of hygiene
- < Full slot-width clip >
  - Delivers exceptional rotation control
- < Reconvertible passive clip design >
  - Provides low friction sliding mechanics of traditional molar tubes while allowing easy wire insertion in the presence of 2nd molars

Page 13-6

## Treatment Versatility

Empower is the first in the industry to offer you the versatility of both interactive and passive bracket designs in one unified system with coordinated in/outs.

Empower's Interactive and Passive bracket designs allow unprecedented treatment options.

### Interactive Bracket

- > Lower friction sliding mechanics early in treatment: exceptional torque and rotation control during working and finishing phases
- > Available in metal upper and lower S-S and in ceramic upper S-S

#### Initial Leveling and Alignment


- Selective Engagement™ clip does not actively engage smaller round and rectangular wires
- Allows smaller wires to unravel debris on for a more comfortable experience and with less clip friction during early treatment

Wires not actively engaged by clip:  
 018 slot: 018 and .016 x .016 and smaller  
 022 slot: 018 and .016 x .022 and smaller

#### Working Stage and Finishing

- Selective Engagement clip actively engages larger wires for precise tooth control during mid-treatment and finishing phases

Wires actively engaged by clip:  
 018 slot: 016 x .022 and larger  
 022 slot: 017 x .025 and larger



Empower's interactive self-ligating system offers the benefits of less friction at the beginning of treatment and more control in the later stages to assist with final torques and rotations.

Empower's passive self-ligating system offers the benefits of less friction and lower force mechanics throughout treatment. No wires are actively engaged by the clip.

> Anterior teeth - interactive brackets  
 low friction, low force during initial leveling and alignment stages,  
 increased torque and rotation control during finishing

> Posterior teeth - passive brackets  
 lower friction mechanics throughout treatment

Matching in/outs between interactive and passive designs means no 1st order wire bending compensations

### Passive Bracket

- > Lower friction sliding mechanics throughout treatment
- > Available in metal upper and lower S-S



## Bond with confidence.

Empower's pads and bases give you the bonding strength you demand and help reduce the risk of emergency appointments due to debonding.




### Quad-Matte™ Base

The patented mechanical lock base of Empower Clear incorporates alumina particles only on the center of the base. This allows for a strong bond in the center of the pad and weaker bond on the edges for easy and predictable debonding.



### Maximum Retention™

The Maximum Retention pad on Empower metal brackets and tubes features photoetched pockets beneath BC gauge mesh to increase pad surface area and ensure safe, secure bonds (see cutaway view).



Empower metal brackets offer multiple textured pad options to meet every need.



**Standard, non-offset pad**  
Slightly larger pad foot allows better fit and grip with steel clips for easy engagement.



**True offset pad**  
Bracket pads offset toward the lingual allows easier bonding on shorter clinical crowns and also works well when the wire bonding.



**Beveled offset pad**  
Slightly larger pad surface area for greater bond strength.

Phone: 1-800-305-5630  
5630-2715  
5630-8825


## Easy open. Easy close.

Significantly faster wire changes than traditionally ligated brackets.<sup>1</sup> The entire Empower family utilizes one double-ended opening instrument providing two opening methods.


<sup>1</sup> Use care not to use the facial notch because over-opening can deform the clips. Use a larger instrument to push the clips inward to prevent the clips from over-opening occlusally.

<sup>1</sup> Chan SS et al. Clin Oral Implants Res. 2010; 21(1): 1-7. Systematic review of self-ligating brackets. AMJ Orthod Dentofacial Orthop. 2010; 137(7): e1-725 e178.


### Original Opening Method™ for Brackets




### Facial Opening Method™ for Brackets



### Facial Opening Method™ for Tubes



### Closing Clip



Close clips easily with simple finger pressure



Page 13-8

## The power to choose.

Empower's breadth of line is unequalled. You get the most features, prescription choices, and treatment possibilities all in one quality, unified system.

More choices. Transforms power.

	Empower™	Empower™ self-aligning	Empower™ 3.0™
Interactive design	•	•	•
Passive design	•	•	•
Maximum Retention™ bonding pads	•	•	•
Quick-Mate™ bonding pads	•	•	•
Standard nonoffset bicuspid pads	•	•	•
True offset bicuspid pads	•	•	•
Extended offset bicuspid pads	•	•	•
Optional hooks available	•	•	•
Individual color coding	•	•	•
Visual Positioning Aids (VPA)	•	•	•
McLaughlin, Bennett, "renew" system	•	•	•
Robt® system	•	•	•
Modified Damon® system	•	•	•
Raincoat® system	•	•	•
Conolly 3.0™s oral vertical slot system	•	•	•
Additional torque options	•	•	•

The American Orthodontics version of the McLaughlin, Bennett, "renew" System is not claimed to be a duplication of any other, nor does American Orthodontics imply that it is endorsed in any way by Dr. McLaughlin, Bennett, or their firm.

The American Orthodontics version of the Robt System is not claimed to be a duplication of any other, nor does American Orthodontics imply that it is endorsed in any way by Dr. Robt.

The American Orthodontics version of the Damon System is not claimed to be a duplication of any other, nor does American Orthodontics imply that it is endorsed in any way by Dr. Damon.

The American Orthodontics version of the Raincoat System is not claimed to be a duplication of any other, nor does American Orthodontics imply that it is endorsed in any way by Dr. Raincoat.



Page 13-9



3524 Washington Avenue Sheboygan, WI USA 53081  
920-457-5051 USA and Canada: 1-800-558-7687 info@americanortho.com www.americanortho.com  
© 2012 American Orthodontics Corporation

**Proudly made in the USA**



LITS1-37 REV 1

Page 13-10





Page 13-11



## CLEAR, BEAUTIFUL BRACES. CLEAR, BEAUTIFUL SMILE.

Presenting Empower Clear Braces – a technologically advanced, ceramic bracket system that will bring you the confidence of a beautiful smile both during and after treatment.

Empower Clear Braces combine two of the orthodontic industry's leading trends - clear cosmetic appliances and self ligation. They are an integral part of Empower, the most complete self ligating bracket system in the world.

The self ligating technology of the Empower system uses integrated clips to hold your arch wire in place, instead of traditional rubber bands (ligatures) that tend to stain and wear down over time. These clips allow your orthodontist to make quicker wire adjustments and deliver the rewarding experience you deserve on the path to your ideal smile.

CLEAR.  
COMFORTABLE.  
PERFORMANCE.

There's no compromise when it comes to how Empower Clear Braces are built. Each bracket is manufactured with leading edge Ceramic Injection Molding technology and aesthetically treated, high performance clips. Empower Clear combines beauty and finesse with performance and precision.

○ Blends beautifully with teeth

Empower Clear Braces are made of strong ceramic material that looks great on your teeth. You'll be noticed, your braces won't.

○ Quicker, easier wire adjustments

State of the art self ligating clips eliminate rubber bands (ligatures) to hold your arch wire in place, allowing your orthodontist to make quicker wire adjustments. This could mean shorter appointments, and may even reduce the total number of visits required during your treatment.

○ Comfortable inside your mouth

Empower Clear's contoured edges and sculpted, low profile design mean a smoother surface against your cheeks and lips.

○ No staining

Self ligating clips eliminate plastic ligatures, which are prone to discoloration and fatigue over time. This, combined with the extremely stable ceramic bracket body material, reduces any risk of staining. Your Empower Clear Braces will stay beautiful throughout treatment.



Page 13-13

## THE PROMISE OF A BETTER SMILE

Your orthodontist is committed to providing you with the brilliant, healthy smile you desire. Empower Clear Braces provide the right tools to make it happen. Leading edge design, aesthetics, and functionality, combined with your orthodontist's clinical expertise help you achieve a smile that will last a lifetime.



Page 13-14

For more information scan this code  
with your mobile phone's QR reader.



AMERICAN  
ORTHODONTICS

[empowerbraces.com/clear](http://empowerbraces.com/clear)

Sheboygan, WI USA 53081

920-457-5051 USA and Canada: 1-800-558-7687

[info@americanortho.com](mailto:info@americanortho.com) [www.americanortho.com](http://www.americanortho.com)

© 2012 American Orthodontics Corporation

**Proudly made in the USA**

LITS1-38

Page 13-15

American Orthodontics 510(k) Submission | 2012

---

**SECTION 14:**  
STERILIZATION / SHELF LIFE

14



**Page 14-1**

**Sterilization / Shelf Life**

**\*\*This section is not applicable. \*\***

American Orthodontics 510(k) Submission | 2012

---

**SECTION 15:**  
**BIOCOMPATIBILITY**





## Biocompatibility

Page 15-1

### Device Information:

Trade Name: Empower Clear  
Common Name: Ceramic Brackets  
Classification Name: Bracket, Ceramic, Orthodontic  
Classification Code: NJM  
Regulation Number: 872.5470

### Equivalent Legally Marketed Devices Information:

Trade Name: Radiance  
Common Name: Ceramic Brackets  
Classification Name: Bracket, Ceramic, Orthodontic  
Classification Code: NJM  
Regulation Number: 872.5470  
510(k) Number: K080749

Trade Name: In-Ovation C  
Common Name: Ceramic Brackets  
Classification Name: Bracket, Ceramic, Orthodontic  
Classification Code: NJM  
Regulation Number: 872.5470  
510(k) Number: K060837

### Summary:

The material, Alumina Oxide [Al<sub>2</sub>O<sub>3</sub>], is used in the manufacturing of the Radiance brackets and In-Ovation C brackets – which is the same material used for Empower Clear. The material was chosen due to the historical use and shown biocompatibility and resistance to corrosion in the oral/human body environment.

Since the Alumina Oxide [Al<sub>2</sub>O<sub>3</sub>] material is chemically inert, is not a metal and is non-conductive, there is no corrosion concerns associated with the bracket relative to the clip or archwire material. Based on these facts, there are no undue risks to patients, users or other persons arising from incompatibilities among the materials used in the manufacture of Empower Clear brackets or among the materials into which the Empower Clear brackets would come into contact with during normal use.

American Orthodontics 510(k) Submission | 2012

---

**SECTION 16:**  
**SOFTWARE**

16



**Page 16-1**

**Software**

**\*\*This section is not applicable. \*\***

American Orthodontics 510(k) Submission | 2012

---

**SECTION 17:**  
**ELECTROMAGNETIC COMPATIBILITY /**  
**ELECTRICAL SAFETY**

17



**Page 17-1**

**Sterilization / Shelf Life**

**\*\*This section is not applicable. \*\***

American Orthodontics 510(k) Submission | 2012

---

**SECTION 18:**  
**PERFORMANCE TESTING – BENCH**



**Page 18-1**

**Performance Testing – Bench**

**\*\*This section is not applicable. \*\***

American Orthodontics 510(k) Submission | 2012

**SECTION 19:**  
**PERFORMANCE TESTING – ANIMAL**





**Page 19-1**

**Performance Testing – Animal**

**\*\*This section is not applicable. \*\***

American Orthodontics 510(k) Submission | 2012

---

**SECTION 20:**  
**PERFORMANCE TESTING – CLINICAL**



**Page 20-1**

**Performance Testing – Clinical**

**\*\*This section is not applicable. \*\***