Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017 K12386Z



MAR 1 2 2013

P.O. Box 708 Warsaw, IN 46581-0708 574 267-6131

13

a

510(k) Summary of Safety and Effectiveness

Sponsor:	Zimmer, Inc. P.O. Box 708 Warsaw, IN 46581-0708		
Contact Person:	Joanna L. Surma Associate Project Manager, Regulatory Affairs Telephone: (574) 371-1642 Fax: (574) 372-4605		
Date:	13 December 2012		
Trade Name:	Zimmer Nexel Total Elbow		
Common Name:	Total Elbow Prosthesis		
Classification	21 CFR § 888.3150 Elbow joint metal/polymer constrained cemented prosthesis		
Product Code:	JDC – Prosthesis, Elbow, Constrained, Cemented		
Predicate Device:	Coonrad/Morrey Total Elbow, manufactured by Zimmer, K001989, cleared 25 July 2000		
	Coonrad/Morrey Total Elbow, manufactured by Zimmer, <u>K053189</u> , cleared 9 December 2005		
Device Description:	The Zimmer Nexel total Elbow is a total elbow prosthesis designed for use with bone cement. It is available in multiple sizes and in right and left configurations.		
	How the Device Functions: The Zimmer Nexel Total Elbow is an implant designed to replace the articulating surfaces of and restore motion to the humeroulnar (elbow) joint. The implant is a constrained device assembly and consists of the following components: humeral component, ulnar component, humeral bearing-A, ulnar bearings-B, axle pin and humeral screws.		
	The Following Scientific Concepts, Design Features and Physical Properties form the Basis for the Zimmer Nexel Total Elbow: The humeral component has a humeral stem designed to be implanted with bone cement into the patient's humeral medullary canal, an anterior flange designed to accept a		

510(k) Summary of Safety and Effectiveness – Zimmer Nexel Total Elbow

bone graft and limit torsional and posterior migration, a humeral yoke with rounded corners to avoid the creation of stress risers within the medial and lateral humeral supracondylar columns, and plasma spray region to enhance fixation to bone cement within the medullary canal, and to improve fatigue strength. The ulnar component has an ulnar stem designed to be implanted with bone cement into the patient's ulnar medullary canal, an ulnar eye that is both highly polished and nitrogen-enriched to limit wear of the apposing polymer bearings, and plasma spray region to enhance fixation to bone cement within the medullary canal. Bearings A and B are designed to broadly distribute joint reaction forces.

Materials Used: The humeral and ulnar components are made of a titanium alloy, the bearings A and B are made of Vitamin E highly cross-linked polyethylene, and the humeral screws and axle pin are made of a cobalt-chromium-molybdenum alloy.

Comparison to the Predicate: The proposed device (Zimmer Nexel Total Elbow) and the predicate device (Coonrad/Morrey Total Elbow) have the same intended use and similar indications for use. The proposed device humeral and ulnar components are very similar in terms of materials used and design/dimensions to the predicate device humeral and ulnar components. The primary differences between the proposed and predicate devices are that the proposed device bearing components-are made of Vitamin E highly cross-linked polyethylene (the predicate device bearing components are made of ultra-high molecular-weight polyethylene), and the proposed device has bearings that articulate on both the outer and inner diameters of the ulnar eye (the predicate device has bearings that articulate on the inner diameter of the ulnar eye). Additionally, the proposed device uses screws and a 1-part axle pin to fix the bearings in place, while the predicate device uses a 2-part (snap-fit) axle pin.

Intended Use:

Indications for use include:

- Elbow joint destruction which significantly compromises the activities of daily living
- Post-traumatic lesions or bone loss contributing to elbow instability
- Ankylosed joints, especially in cases of bilateral ankylosis from causes other than active sepsis
- Advanced rheumatoid, post-traumatic, or degenerative arthritis with incapacitating pain
- Instability or loss of motion when the degree of joint or soft tissue damage precludes reliable osteosynthesis
- Acute comminuted articular fracture of the elbow joint surfaces that precludes less radical procedures, including 13-C3 fractures of the distal humerus
- Revision arthroplasty

Caution: This device is intended for cemented use only.

510(k) Summary of Safety and Effectiveness – Zimmer Nexel Total Elbow

Comparison to Predicate Device:

Performance Data:

The Zimmer Nexel Total Elbow is substantially equivalent to the predicate device in terms of form and function. The Zimmer Nexel Total Elbow and the predicate device share similar intended uses and indications for use.

Non-Clinical Performance Testing Conducted:

- Stem Fatigue Testing
- Wear Testing
- Durability Testing
- Modular Connection Fatigue Testing

Non-Clinical Performance Testing Conclusions:

Non-clinical testing demonstrated that the New Zimmer Total Elbow meets performance requirements as defined by Design Control activities and is substantially equivalent to the predicate device in terms of safety and efficacy.

In this case, clinical data and conclusions were not needed to demonstrate substantial equivalence.

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Letter dated: March 12, 2013

Food and Drug Administration 10903 New Hampshire Avenue Document Control Center – WO66-G609 Silver Spring, MD 20993-0002

Zimmer, Incorporated % Ms. Joanna L. Surma Associate Project Manager, Regulatory Affairs P.O. Box 708 Warsaw, Indiana 46581-0708

Re: K123862

Trade/Device Name: Zimmer Nexel Total Elbow Regulation Number: 21 CFR 888.3150 Regulation Name: Elbow joint metal/polymer constrained cemented prosthesis Regulatory Class: II Product Code: JDC Dated: December 13, 2012 Received: December 14, 2012

Dear Ms. Surma:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

Page 2 – Ms. Joanna L. Surma

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <u>http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm</u> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm,

Sincerely yours,

Mark N.Melkerson -S

Mark N. Melkerson Director Division of Orthopedic Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K123862

Device Name:

Zimmer Nexel Total Elbow

Indications for Use:

Indications for use include:

- Elbow joint destruction which significantly compromises the activities of daily living
- Post-traumatic lesions or bone loss contributing to elbow instability
- Ankylosed joints, especially in cases of bilateral ankylosis from causes other than active sepsis
- Advanced rheumatoid, post-traumatic, or degenerative arthritis with incapacitating pain
- Instability or loss of motion when the degree of joint or soft tissue damage precludes reliable osteosynthesis
- Acute comminuted articular fracture of the elbow joint surfaces that precludes less radical procedures, including 13-C3 fractures of the distal humerus
- Revision arthroplasty

CAUTION: This device is intended for cemented use only.

Prescription Use <u>X</u> (Part 21 CFR 801 Subpart D) AND/OR

Over-The-Counter Use (21 CFR 807 Subpart C)

(Please do not write below this line - Continue on another page if needed)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Casey L. Hanley, Ph.D.

Page 1 of 1.

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration 10903 New Hampshire Avenue Document Control Center – WO66-G609 Silver Spring, MD 20993-0002

Zimmer, Incorporated % Ms. Joanna L. Surma Associate Project Manager, Regulatory Affairs P.O. Box 708 Warsaw, Indiana 46581-0708 Letter dated: March 12, 2013

Re: K123862

Trade/Device Name: Zimmer Nexel Total Elbow Regulation Number: 21 CFR 888.3150 Regulation Name: Elbow joint metal/polymer constrained cemented prosthesis Regulatory Class: II Product Code: JDC Dated: December 13, 2012 Received: December 14, 2012

Dear Ms. Surma:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

Page 2 - Ms. Joanna L. Surma

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <u>http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm</u> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <u>http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm</u>.

. Sincerely yours,

Mark N. Melkerson -S

Mark N. Melkerson Director Division of Orthopedic Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

Page 3 – Ms. Joanna L. Surma

Concurrence & Template History Page [THIS PAGE IS INCLUDED IN IMAGE COPY ONLY]

Full Submission Number: K123862

For Office of Compliance Contact Information:

http://insideportlets.fda.gov:9010/portal/page?_pageid=197,415881&_dad=portal&_schema=PORTAL&org=318

For Office of Surveillance and Biometrics Contact Information:

http://insideportlets.fda.gov:9010/portal/page?_pageid=197,415881&_dad=portal&_schema=PORTAL&org=423

Digital	Signature Concurrence Table
Reviewer Sign-Off	Peter Allen 3/10/13
Branch Chief Sign-Off	Casey Hanley
Division Sign-Off	Mark N. Melkerson -S 2013.03.1217:14:06-04'00'

f/t:PGA:tmj:3/12/13:eaf4:3/12/13

Template Name: K1(A) - SE after 1996

Template History:

Date of Update	By	Description of Update
7/27/09	Brandi Stuart	Added Updates to Boiler Table
8/7/09	Brandi Stuart	Updated HFZ Table
1/11/10	Diane Garcia	Liability/Warranty sentence added at bottom of 1 st page
10/4/11	M. McCabe Janicki	Removed IFU sheet and placed in Forms
9/25/12	Edwena Jones	Added digital signature format
12/12/12	M. McCabe Janicki	Added an extra line between letter signature block and the word
		"Enclosure". Also, added a missing digit in 4-digit extension on
		letterhead zip code: "002" should be "0002".

Indications for Use

510(k) Number (if known): <u>K123862</u>

Device Name:

Zimmer Nexel Total Elbow

Indications for Use:

Indications for use include:

- Elbow joint destruction which significantly compromises the activities of daily living
- Post-traumatic lesions or bone loss contributing to elbow instability
- Ankylosed joints, especially in cases of bilateral ankylosis from causes other than active sepsis
- Advanced rheumatoid, post-traumatic, or degenerative arthritis with incapacitating pain
- Instability or loss of motion when the degree of joint or soft tissue damage precludes reliable osteosynthesis
- Acute comminuted articular fracture of the elbow joint surfaces that precludes less radical procedures, including 13-C3 fractures of the distal humerus
- Revision arthroplasty

CAUTION: This device is intended for cemented use only.

Prescription Use X (Part 21 CFR 801 Subpart D) AND/OR

Over-The-Counter Use (21 CFR 807 Subpart C)

(Please do not write below this line - Continue on another page if needed)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Casey L. Hanley-Ph.D. Division of Orthopaedic Devices

Page 1 of 1

REASON FOR ERROR E-11 HANG UP OR LINE

E-2) BUSY E-4) NO FACSIMILE CONNECTION



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration 10903 New Hampshire Avenue Document Control Center - WO66-G609 Silver Spring, MD 20993-0002

Zimmer, Incorporated % Ms. Joanna L. Surma Associate Project Manager, Regulatory Affairs P.O. Box 708

Letter dated: March 12, 2013

Warsaw, Indiana 46581-0708

Re: K123862

Trade/Device Name: Zimmer Nexel Total Elbow Regulation Number: 21 CFR 888.3150 Regulation Name: Elbow joint metal/polymer constrained comented prosthesis Regulatory Class: 11 Product Code: JDC Dated: December 13, 2012 Received: December 14, 2012

Dear Ms. Surma:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017

		COVER S	HEET	MEMORAI	N D U M	Food and D Office of D Office of I	evice Eva	aluation &
From:	Rev	viewer Name	Peter /	Allen	<u> </u>			
Subject:	510	(k) Number	K1238	62			_	
То:	The	Record						
o Refuse http://er 202%20 o Hold (/	ed to com.1 007.de Additi	ida.gov/eRoomRec oc) onal Information	is conside /Files/CDRH or Telephor	<u>H3/CDRHPremarke</u> ne Hold).	w cycle, See Screening Cho t <u>Notification510kProgram/0_5</u> e below), Withdrawn, etc.).		ing%20C	hecklist%2079
	No	t Substantially Ec	quivalent (N	ISE) Codes				
	0000 00000	NO NI NQ NU NP NS NL NM NC NH	NSE for ne NSE for ne effective NSE for lac NSE no re NSE for lac NSE pre-a NSE post-a	ew intended use A iness ck of performance sponse ck of performance mendment device	e data AND no response e call for PMAs (515i) ce requires PMAs			
	0	TR	NSE for tra	ansitional device				
Post St. Mill Post - Post	-		or a final cle	earance decision	(i.e., SE, SE with Limitation	s, etc.):	YES	NO
		Use Page			Attach IFU		X	
		ry /510(k) Statem	· · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	Attach Summary	I Decision	X	
AL CONTRACTORISTIC TROUBLE CLASSE		curate Statemen	IL.		Must be present for a Fina	i Decision	х	
Is the dev		n include Class III	l Summon	<u>^</u>	Must be present for a Fina	Decision		x
Does firm	refer s, plea	ence standards?	en e	NARE OF A CONTRACT OF A SUB-	om/morechoices/fdaforms/	···· ;	x	
(Pleas http://e	e spe room	nation product? ecify category .fda.gov/eRoomRe %20PRODUCT%2	, se q/Files/CDR 20ALGORITH	e H3/CDRHPremarke HM%20(REVISED%	etNotification510kProgram/0_ 6203-12-03).DQC	413b/CO		x
(Guida Repro	ance cessi	ed Single-Úse Me	DA Staff – edical Devic	ces, <u>http://www.fd</u>	lation Data in 510(k)s for a.gov/cdrh/ode/guidance/12	2 <u>16.html</u>)		X
	A. 14	ntended for pedia				» i istos toto		X .
Is this a p	rescr	iption device? (If	ooth prescr	PM EDA 2674 C	eck both boxes.) ertification with Requireme	nte of	X	1 <u>2000</u> 1 1
		ov Data Bank?						X

For United States-based clinical studies only: Did the application include a completed FORM FDA 3674, Certification with Requirements of ClinicalTrials.gov Data Bank? (If study was

Is clinical data necessary to support the review of this 510(k)?


Rev. 9/20/12 - added digital concurrence table

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

х

21 CFR 888.3150	(*If unclassified, see 510	JDC	
Regulation Number	Class*	Product Code	
Is this device subject to the Tra Guidance, <u>http://www.fda.g</u>	acking Regulation? (Medical I ov/cdrh/comp/guidance/169.h		C. X
Companion Diagnostic			X
Device Contains Battery			x
MR Conditional			x
Mobile Application			x
Nanotechnology		-	x
Transitional Adolescent B (18 - old)	<= 21; No special considerati	ons compared to adults => 21 yea	ars X
group, different from adults age procedures, etc.)	e ≥ 21 (different device desig	· ·	x
Adolescent (12 years -< 18 yea	irs old)		X
Child (2 years -< 12 years old)			· x
Infant (29 days -< 2 years old)			x
Neonate/Newborn (Birth to 28)	days)		x
All Pediatric Patients age<=21			x
Does this device include an An	imal Tissue Source?		x
conducted in the United States applicant must be contacted to			i i

Digital	Signature Concurrence Table
Peter G. Allen 2013.03.11 09:13:11 -04	
Branch Chief Sign-Off	Casey Hanley - S 2013.03.12 09:22:09 -04'00'
Division Sign-Off	Mark N. Melkerson -S 2013.03.12 17:12:09 -04'00'



* 510(k) Submissions compare new device to marketed devices. FDA requests additional Information if the relationship between marketed and "predicate" (pre-Amendments or reclassified post-Amendments) devices is unclear.

** This decision is normally based on descriptive information alone, but limited testing information is sometimes required.

*** Data may be in the 510(k), other 510(k)s, the Center's classification files, or literature.



DEPARTMENT OF HEALTH AND HUMAN SERVICES

MEMORANDUM

Food and Drug Administration Office of Device Evaluation 10903 New Hampshire Avenue Silver Spring, MD 20993

Premarket Notification [510(k)] Review Traditional/Abbreviated

<u>K123862</u>

Date: March 8, 2013Office: ODETo: The RecordDivision: DODFrom: Peter G. Allen, M.S.; Biomedical EngineerBranch: JFOB510(k) Holder: Zimmer, Inc.Device Name: Zimmer Nexel Total ElbowContact: Ms. Joanna Surma, MS, RAC, Regulatory Affairs, Associate Project ManagerPhone: 574-371-1642Fax: 574-372-4605Email: joanna.surma@zimmer.com

I. Purpose and Submission Summary

The 510(k) holder would like to introduce Zimmer Nexel Total Elbow into interstate commerce. The Nexel Total Elbow is a constrained elbow prosthesis intended for cemented use. This is a new device to the U.S. market.

The sponsor provided, via interactive review, additional information including EtO residual testing results, details on colorants used in their instrumentation, updated labeling to include appropriate MR statements, additional discussion regarding the differences in range of motion between their device and the predicate device, edits to all relevant parts of the submission to replace the device descriptive term "semi-constrained" with the appropriate term "constrained", and revisions to the 510(k) Summary.

The materials, design, sizes, indications, intended use, sterilization, and packaging, for the subject device are similar and/or identical to legally marketed predicate devices. The sponsor has provided adequate characterization in terms of device design, size, material, intended used and testing to demonstrate that the subject device components are equivalent to marketed predicate devices. <u>Therefore, I recommend SE</u>.

II. Administrative Requirements

	Yes	No N
Indications for Use page (Indicate if: Prescription or OTC)	x	
Truthful and Accuracy Statement	x	
510(k) Summary or 510(k) Statement	x	
Standards Form	×	


III. Device Description

	Yes	No	N/A
Is the device life-supporting or life sustaining?		x	

1

	Yes,	Nö	N/A
Is the device an implant (implanted longer than 30 days)?	x		ļ
Does the device design use software?		X	
Is the device sterile?	x		
Is the device reusable (not reprocessed single use)?		×	
Are "cleaning" instructions included for the end user?		1	×

The Zimmer Nexel Total Elbow is designed to replace the articulating surfaces of and restore motion to the humeroulnar (elbow) joint. The implant is a constrained device assembly and consists of a humeral component, an ulnar component, a humeral bearing-A, 2 ulnar bearings-B, an axle pin and 2 humeral screws.



Exploded view

Assembled view

Humeral Component

The humeral component is made from wrought *Tivanium* alloy (Ti-6AI-4V, ASTM F136) and has the following major design features: a stem, an anterior flange, a yoke and plasma spray regions.

The humeral stem is intended to be implanted with bone cement into the patient's humeral medullary canal. It has a triangular cross-section, which is designed to prevent rotation of the humeral component within the patient's medullary canal; and, it is straight along its proximal-distal axis to conform to the straight shaft of the distal humeral medullary canal.

The humeral anterior flange is designed to accept a bone graft, which is wedged between the anterior flange and the exterior cortical bone of the humeral shaft. This anterior flange and bone graft combination are intended to improve both antirotational stability and posterior migration of the humeral component, thereby resisting stem loosening.

The humeral yoke has rounded corners (as opposed to angular) on its proximal surface in order to prevent stress risers (that are often associated with bone preparation cuts required for angled corners abutting bone). This rounded yoke design is intended to preserve the strength of (and avoid fracture of) the medial and lateral supracondylar columns. The yoke walls also house partially threaded holes with Spiralock technology, which

engage with the humeral screws.

The humeral component has several regions covered in plasma spray (which is contained within the Plasma Sprayed Titanium-6Ai-4V Alloy Coating Material Master File, FDA Ref: MAF-1909, Letter of Authorization can be found in Section 21). The plasma spray regions are designed to increase the coefficient of friction between the implant and the bone cement, thereby improving anti-rotational stability and aiding in resistance to movement of the implant with respect to the humerus.

In addition to the aforementioned major design features, the humeral component also includes the following features to accept the axle-pin and bearings A and B: v-grooves at the base of the keyways to accept the axle-pin, a hole to accept humeral bearing-A and two keyways to accept ulnar bearings-B.

The Nexel Total Elbow was designed to provide at least as much of the necessary flexion and v/v laxity as required of a normal, healthy elbow while maximizing the thickness and optimizing the location and form of the Vitamin-E UHMWPE Bearings to provide improved wear resistance and longer-term load-bearing capability as compared to the predicate device. Thus, any flexion and v/v angulation deemed unnecessary for a full healthy elbow range of motion was traded for optimization of device performance; this is further substantiated by the results of the head-to-head Wear and Durability tests, wherein both (Nexel and predicate) devices were subjected to high-loading over large flexion arcs with induced v/v angulation to simulate clinically-relevant, yet aggressive, long-term use in-vivo.

Ulnar Component

The ulnar component is made from wrought *Tivanium* alloy (Ti-6AI-4V, ASTM F136) and has the following major design features: a stem, an eye, and a plasma spray region.

The ulnar stem is intended to be implanted with bone cement into the patient's ulnar medullary canal. It has a quadrangular cross-section, which is designed to prevent rotation of the ulnar component within the patient's medullary canal; and, it is curved along its proximal-distal axis to conform to the curved shaft of the proximal ulnar medullary canal.

The ulnar eye is both highly polished and nitrogen-enriched (as recommended in Section J, Titanium Articular Surface Coatings, of *Class II Special Controls Guidance Document: Knee Joint Patellofemorotibial and Femorotibial Metal/Polymer Porous-Coated Uncemented; Guidance for Industry and FDA)* to limit wear of the opposing Bearings A and B upon which it articulates. This a nitrogen ion implantation (also known as nitriding) process to increase titanium alloy surface hardness on the region specified in the engineering drawings (both inner and outer articulating surfaces of the ulnar eye – as clarified in phone call with sponsor on 2/22/13 and their email of 2/27/13). The articulating surfaces of the ulnar component are polished to an Ra value of (b)(4) or less and the articulating surfaces of the bearings A and B are machined to an Ra value of (b)(4) or less to provide for smooth articulating surfaces.

The proximal portion of the ulnar stem is covered with plasma spray (which is contained within the Coatings Characterization of Orthopedic Titanium 6A/4V Coatings Material Master File, FDA Ref: MAF-652, Letter of Authorization can be found in Section 22). The plasma spray region is designed to increase the coefficient of friction between the implant and the bone cement, thereby improving anti-rotational stability and aiding in resistance to movement of the ulnar stem with respect to the ulna.

The ulnar component is side-specific, meaning that there is a left ulnar component and a right ulnar component available for every size offering.

Humeral Bearings-A

The Humeral Bearing-A is made from Vivacit-E (Vitamin E [α -tocopherol] stabilized, highly crosslinked ultra-high molecular weight polyethylene), which is contained within the Vitamin E Stabilized HXPE Material Master File authored by Zimmer (FDA Ref. MAF-1868). Humeral bearing-A is designed to articulate against the ulnar

component throughout most of the humeroulnar joint range of motion, and has the following major design features: concave articulation surface, anterior and posterior ridges, a fixation peg, four peripheral bumps, and medial and lateral rails.

The humeral bearing-A concave articulation surface is intended to, along with the articulation surface of the ulnar bearings-B, allow for load sharing from the ulnar component onto humeral bearing-A throughout most of the range of motion of the humeroulnar joint. The posterior ridge specifically supports high-flexion motions.

The humeral bearing-A fixation peg is designed to mate to the hole at the proximal base of the humeral yoke with a press (not a snap) fit, which is intended to limit motion between humeral bearing-A and the humeral yoke.

The humeral bearing-A has 4 peripheral bumps that are designed to mate against the interior walls of the humeral yoke with a press-fit, which is intended to limit motion (and therefore backside wear) of the humeral bearing-A relative to the humeral yoke.

The humeral bearing-A medial and lateral rails are designed to allow clearance for (as opposed to a direct mate with) ulnar bearings-B, which directly mate to the humeral yoke and to the axle pin.

Ulnar Bearings-B

The ulnar bearings-B (two per implant assembly) are made from Vivacit-E (Vitamin E [α-tocopherol] stabilized, highly crosslinked ultra-high molecular weight polyethylene), which is contained within the Zimmer Vitamin E Stabilized HXPE Material Master File (FDA Ref: MAF-1868). Ulnar bearings-B are designed to articulate against the ulnar component throughout the entire humeroulnar joint range of motion, and have the following major design features: articulation surface, anti-rotation key, through-hole and compressible semi-circle.

The ulnar bearings-B articulation surface is intended to, along with the articulation surface of the ulnar bearing-A, allow for load sharing from the ulnar component onto humeral bearing-B throughout the range of motion of the humeroulnar joint.

The ulnar bearings-B anti-rotation keys mate into the humeral yoke keyways with a press-fit to limit translational and rotational motion of ulnar bearings-B with respect to the humeral component.

The compressible semi-circle of one ulnar bearing-B is designed to align with the compressible semi-circle of the second ulnar bearing-B such that a nearly full circle is created. This design feature is intended to ensure proper alignment of ulnar bearings-B with respect to each other (during assembly of ulnar bearings-B to the axle pin using the ulnar bearing assembly tool) and, after proper assembly to the axle pin, to the humeral yoke keyways. Due to the fact that the compressible semi-circle is more easily deformed than the broad, apposing faces of ulnar bearings-B, it creates a press-fit between the ulnar bearings-B and between each ulnar bearing-B and the opposing interior walls of the humeral yoke. This press-fit is intended to limit mediolateral motion of the ulnar bearings-B along the axle-pin, thereby limiting backside wear of the ulnar bearings-B.

The ulnar bearings-B through-hole mates around the axle pin with a press-fit intended to limit motion of ulnar bearings-B with respect to the axle pin, and thus limit backside wear of ulnar bearings-B. The position of the through-hole is eccentric (with respect to the outer diameter of ulnar bearings-B) in order to place the largest cross-section of Vivacit-E material in opposition to the highest loads (joint reaction forces) induced throughout the flexion-extension cycle of the humeroulnar joint.

Axle Pin

The axle pin is made from Zimaloy (Cobalt-Chromium-Molybdenum Alloy, ASTM F75), is designed to remain in a fixed position (not translate nor rotate) throughout the entire humeroulnar joint range of motion, and has the following major design feature: v-groove.

The v-groove is designed to mate with the taper of the humeral screws that, upon screw thread engagement in the humeral yoke, rigidly compresses the axle pin against the humeral yoke. Additionally, the axle pin is designed to rigidly mate with ulnar bearings-B, the humeral screws and the humeral yoke, creating a rigid assembly along which joint reaction forces are intended to be shared.

Humeral Screws


The humeral screws (two per implant assembly) are made from Zimaloy (Cobalt-Chromium-Molybdenum Alloy, ASTM F75), are designed to remain in a fixed position throughout the entire humeroulnar joint range of motion, and have the following major design features: headless, threads, pilot and taper.

The humeral screw threads engage with threads of the humeral yoke while the humeral screw pilot engages a hole in the wall of the humeral yoke and the humeral screw taper engages with the axle pin v-groove, compressing the ends of the axle pin against the v-groove at the base of the humeral yoke keyway.

The humeral screws are headless, a feature that allows the humeral screws to countersink into the humeral yoke, which is intended to ensure the constant application of a clamp-load to the axle-pin and to avoid soft tissue impingement or irritation associated with a protruding screw head.

While the humeral screws are designed to apply a consistent clamp-load to the axle pin, they are not directly loaded by any compressive joint reaction force during the flexion-extension cycle (this design feature is intended to avoid screw back-out by not allowing any complete transverse load-reversal across the humeral screw threads). Instead, compressive forces are transmitted from the ulnar implant component to the ulnar bearings-B to the axle pin to the y-prooves at the base of the humeral voke keyways and finally throughout the humeral component.

Other



poly bearings (humeral bearing-A and ulgar bearings-B) are machined from (10)(4)

crewe and ayle nine are machin

Zimmer Nexel Total Elbow components are not compatible with any other devices outside the subject Zimmer Nexel Total Elbow system.

All Zimmer Nexel Total Elbow humeral component sizes are designed to mate with all Zimmer Nexel Total Elbow ulnar component sizes. The axle pin size 4, humeral bearing-A size 4, and ulnar bearings-B size 4 are designed to mate with all size 4 humeral components. The axle pin size 5/6, humeral bearing-A size 5/6, and ulnar bearings-B size 5/6 are designed to mate with all size 5 and size 6 humeral components. The humeral screws are designed to mate with all Zimmer Nexel Total Elbow humeral components and axle pin sizes.

Representative engineering drawings are provided in the Device Description section of the submission.

The sponsor has provided a table of class I and class II instrumentation applicable to the use of this device. Materials include predominantly stainless steel, as well as Radel, Ultern, LDPE, and silicone.

Two instruments were identified as having colorant additives in them. These were an Ultern screw driver handle and a silicone screw removal driver handle. The screw driver is a single-use only disposable device. The same colorants used in the polymeric instrument components have been used in other Zimmer device systems including K040593 and K933785.

Reviewer Comment: No information on colorants was provided in the original submission. In an interactive review teleconference (and email) with the sponsor on 2/22/13 the sponsor was asked to either identify predicate device systems that utilize these colorants (if any) or provide biocompatibility information. In their email response dated and received on 2/27/13 the sponsor identified that two colorants were used with the subject instrumentation. They also identified two predicate device systems that have utilized these same colorants in similar instrument components. The colorants were not identified by name but since they are just used in instrument handles and have been used for 9 and 20 years, respectively (and one is disposable, single-use only), no additional information will be asked for at this time. Therefore, the response is acceptable; and no

further information is required to 'grandfather' in these colorants. The sponsor was notified, however, that in the future any new colorants will need biocompatibility information provided in the 510(k) of the associated device system.

This device description is now adequate as information on colorants used in the instrumentation has been provided. The device description is adequate to demonstrate the similarities and differences in device design to the referenced predicate devices, and in helping to establish an SE determination based on these comparisons.

IV. Indications for Use

Indications for use include:

- Elbow joint destruction which significantly compromises the activities of daily living.
- Post-traumatic lesions or bone loss contributing to elbow instability.
- Ankylosed joints, especially in cases of bilateral ankylosis from causes other than active sepsis.
- Advanced rheumatoid, post-traumatic, or degenerative arthritis with incapacitating pain.
- Instability or loss of motion when the degree of joint or soft tissue damage precludes reliable osteosynthesis.
- Acute comminuted articular fracture of the elbow joint surfaces that precludes less radical procedures, including 13-C3 fractures of the distal humerus.
- Revision arthroplasty.

CAUTION: This device is intended for cemented use only.

Reviewer Comment: The indications are the same or similar to those identified in the predicate device submissions. The indications are consistent across the Package Insert, 510(k) Summary, and Indications for Use form. This is acceptable.

V. Predicate Device Comparison

K001989, Coonrad/Morrey Total Elbow, Zimmer, Inc. K053189, Coonrad/Morrey Total Elbow, Zimmer, Inc.

Reviewer Comment: The sponsor has identified acceptable predicate devices.

VI. Labeling

Labeling includes draft labels, a draft package insert (IFU), separate IFU for cleaning and sterilization of nonsterile instrumentation (provided as a response to RTA deficiencies), and a draft surgical technique.

Symbols on the outer package labels are paired with appropriate descriptive text. Package labels appear to contain all necessary information including the prescription use only symbol (Rx symbol), sterility, sterilization method, single use only statement, cemented use only statement, caution statement (see instructions for use), material, quantity, expiration date (10 years and 5 years), size, name, lot number, catalogue number, sponsor name, and address. The outer package labels appear appropriate.

A specific package insert for the Nexel Total Elbow was provided which covers the components of this subject submission. It contains all appropriate sections, including sponsor name and contact information, device name, cemented use only statement, prescription use statement, device description, materials, indications, contraindications, warnings, precautions, adverse events, and sterility. Appropriate MR language noting that the Nexel Total Elbow has not been evaluated for safety and compatibility, nor tested for heating or migration in the MR environment, *is* included in the revised package insert, provided via email on 2/27/13. Resterilization of UHMWPE implants is specifically not recommended.

Sterilization/resterilization parameters for metal implant components and instruments are provided in the draft package insert and IFU for cleaning and sterilization of non-sterile instrumentation, respectively. The parameters are in accordance with tables 4 and/or 5 of ANSI/AAMI ST79. See Sterilization section below for details.

The following parameters are recommended for resterilization of metal implants: steam heat using a gravity

displacement cycle at 121° C (250° F) for 30 minutes (minimum), allowing 15 minutes (minimum) of drying time, and steam heat using a gravity displacement cycle at 132° C (270° F) for 15 minutes (minimum), allowing 15 minutes (minimum) of drying time. For instruments a pre-vacuum cycle at 132° C (270° F) for 4 minutes (minimum), allowing 30 minutes (minimum) of drying time, is recommended. These are the same as for other Zimmer systems.

The surgical technique appears acceptable. It provides all detailed aspects of the surgery with device specific instrumentation. It includes schematics and photos in conjunction with the text.

Reviewer Comment: In their 2/27/13 email response to our interactive review deficiencies discussed via phone (and email) on 2/22/13 the sponsor provided an updated package insert that included the above referenced MR statements. The labeling is acceptable.

A revised 510(k) Summary was provided via email on 2/28/13. It includes requested revisions such as a complete device description section (as per the RTA checklist and 510(k) Summary checklist (per 21 CFR 807.92]), corrections to the product code and regulation/classification names, and replacement of the word "semi-constrained" with "constrained". The 510(k) Summary is now acceptable.

In addition, in the 2/27/12 email the sponsor revised the following sections of the submission by replacing the term "semi-constrained" with "constrained", as requested:

- original K123862 submission page 16 see page 1 of Attachment 2
- original K123862 submission page 56 see page 1 of Attachment 4 -
- original K123862 submission page 61 see page 1 of Attachment 5
- original K123862 submission page 74 see page 14 of Attachment 5
- -RTA Response page 4 – see page 2 of Attachment 6
- RTA Response page 27 see page 9 of Attachment 7
- original K123862 submission page 57 see page 2 of Attachment 4

VII. Sterilization/Shelf Life/Reuse

All implant components are provided sterile. The metallic humeral component, ulnar component, and humeral screws are sterilized by gamma irradiation. The articulation kit (humeral bearing-A, 2 ulnar bearings -B, and an axle pin) which contains highly-crosslinked Vitamin E UHMWPE components is sterilized by ethylene oxide. Some instruments are also provided sterile.

Metallic implant components provided sterile by radiation:

Method: gamma radiation from a Co 60 source

Dosage: 25 - 37 kGy (2.5 - 3.7 Mrads) Sterility Assurance Level: 10⁻⁶

Validation Method: ANSI/AAMI/ISO 11137-1; 2006 (Sterilization of health care products. Part 1, Radiation. Requirements for development, validation and routine control of a sterilization process for medical Devices.), 11137-2: 2006 (Sterilization of health care products - Radiation - Part 2: Establishing the sterilization dose), and AAMI TIR 33: 2005, (Sterilization of health care products - Radiation Substantiation of a selected sterilization dose - Method VDmax).

Pyrogenicity: Components will not be labeled as nonpyrogenic or pyrogen free. Per USP XXIII (161), requirements for specified endotoxin levels do not apply to orthopedic implants.

Articulation kit (1 Vivacit-E humeral bearing-A and 2 Vivacit-E ulnar bearings-B and an axle-pin) provided sterile by EO

Method: Ethylene Oxide

Validation: AAMI/ISO 11135-1: 2007 (Sterilization of health care products - Ethylene Oxide Part 1: Requirements for development, validation, and routine control of a sterilization process for medical devices). SAL: 10-6

Pyrogenicity: Components will not be labeled as pyrogen free.

Residual Levels: The maximum levels of residuals of EO and ethylene chlorohydrin that remain on the device do not exceed the maximum allowable limits described in AAMI/ISO 10993-7: 2008 (Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals) for permanent contact devices. The maximum allowable levels are as follows:

- The average daily dose of EO to patient shall not exceed 0.1mg/d. In addition the maximum EO dose shall not exceed 4mg in the first 24 hours; 60mg in the first 30 days; and 2.5 grams in a lifetime.
- The average daily dose of ECH to patient shall not exceed 0.4mg/d. In addition the maximum ECH dose shall not exceed 9mg in the first 24 hours; 60mg in the first 30 days; and 10 grams in a lifetime.

The sponsor provided EO residuals testing results via an interactive review email sent/received on 2/14/13. The results are below the maximum allowable limits noted above for each specified time point.

Class II accessories provided sterile

The following accessories listed in Section 11a (Device Description) of the submission are provided sterile and are intended to be used during only one surgery (single use) and are not intended to be end-user sterilized.

Catalog Numb e r	Description	Sterilization Method	Radiation Dose	Validation Method	SAL
00-8401-088-00	Humeral Screw Installation Driver Bit	Gamma Irradiation	(b)(4)	AAMI TIR33, VD _{MAX}	10 ⁻⁶ or better
00-8401-080-00	Torque-Limiting Humeral Screw Driver	Gamma Irradiation		AAMI 11137-2, VD _{MAX}	10 ⁻⁶ or better

Table 1. Class II Accessories to the Implant provided sterile







Humeral Screw Kit -(b)(4)

2



The above listed materials are identical to those currently in use for Zimmer sterile implant packages with the only exception being that the configuration of the interior components correspond to the particular product being packaged.

The Class II Accessories to the Implant (listed in Table 1, above) that are provided sterile are packaged as follows:

Humeral Screw Installation Driver Bit -

Torque-Limiting Humeral Screw Driver -

(b)(4

The above listed materials are identical to those currently in use for other similar Zimmer sterile implant packages with the only exception being that the configuration of the interior components correspond to the particular product being packaged.

<u>Recommended and validated sterilization parameters for end users of reusable class II accessories</u> (instruments) provided non-sterile:

Method: moist heat (steam)

Validation: AAMI/ANSI/ISO 17665-1:2006 Sterilization of Health Care Products - Moist Heat - Part 1: Requirements for the development, validation, and routine control of a sterilization process for medical devices. Half cycle method (Overkill method).

Sterility Assurance Level: 10⁻⁶

Cycle: dynamic air removal (pre-vacuum) for 4 minutes at 132° C with 30 minutes minimum dry time, using recommended sterilization wraps.

Parameters conform to those in Table 5 of AAMI ST79: 2010.

Shelf-Life:

The shelf life of the packaging for the Zimmer Nexel Total Elbow humeral component, ulnar component and humeral screws (note: the "humeral screws" are also referred to as the "humeral screw kit") is 10 years. The shelf life of the articulation kit (the humeral bearing-A, humeral bearings-B and the axle pin) is 5 years. This testing was performed according to ASTM F1980 (*Standard Guide for Accelerated Aging for Sterile Medical Device Packages*) and ANSI/AAMI/ISO 11607: 2006 (*Packaging for Terminally Sterilized Medical Devices*). The packaging materials demonstrate no degradation over time and the shelf life test data assures sterile packaging integrity.

Reviewer Comment: Sterilization/packaging is acceptable as currently presented, and is identical to earlier predicate submissions of theirs. Resterilization of UHMWPE device implants is specifically not recommended.

VIII. Blocompatibility

The device materials, CoCrMo alloy and titanium alloy, have a long (~40 years) and successful history in orthopedic total joint replacement. The materials also conform to recognized material standards. According to

FDA reviewer Michael Kasser, Ph.D., the Vivacit-E UHMWPE (Vitamin E [o-tocopherol] stabilized, highly crosslinked ultra-high molecular weight polyethylene), which is contained within the Zimmer Vitamin E Stabilized HXPE Material Master File (FDA Ref: MAF-1868) has been previously cleared for use in other Zimmer devices and does not present a biocompatibility concern.

On page 15 of the sponsor's 1/17/13 RTA response it is noted that all of the class II accessories (instruments) are external communicating devices contacting tissue/bone with limited exposure (less than 24 hour duration), according to definitions presented in Section 5 of ISO 10993-1.

No information on colorants was provided in the original submission. In an interactive review teleconference (and email) with the sponsor on 2/22/13 the sponsor was asked to either identify any colorants (if any) and the predicate device systems that utilize these colorants, or provide biocompatibility information. In their email response dated and received on 2/27/13 the sponsor identified that two colorants were used with the subject instrumentation. They also identified two predicate device systems that have utilized these same colorants in similar instrument components (K040593, K933785). The colorants were not identified by name but since they are just used in instrument handles and have been used for 9 and 20 years, respectively (and one is disposable, single-use only), no additional information will be asked for at this time. Therefore, the response is acceptable; and no further information is required to 'grandfather' in these colorants. The sponsor was notified, however, that in the future any new colorants will need biocompatibility information provided in the 510(k) of the associated device system.

Reviewer Comment: No additional biocompatibility information is needed.

IX. <u>Electromagnetic Compatibility and Electrical, Mechanical and Thermal Safety</u> None.

Reviewer Comment: Electromagnetic compatibility data is typically not provided in orthopedic 510(k) submissions. Electrical and thermal safety issues are not relevant (typically) to orthopedic devices. The sponsor's package insert, as amended in the sponsor's 2/27/13 email response, includes appropriate MR language. Mechanical safety issues are addressed in the following section on performance testing.

X. Performance Testing - Bench

The Nexel Total Elbow has been evaluated in comparison to its predecessor and predicate device system, the Coonrad/Morrey Total Elbow. The following table identifies the types of testing conducted by the sponsor. Test results demonstrate that the subject device performed as well as or better than the predicate device for each test. This is as expected due to the very similar nature of the device designs. The subject device incorporates iterative design changes from the predicate, implementing improvements in some areas (e.g., humeral stern fatigue strength, wear, durability) and maintaining equivalence in others (ulnar stern fatigue strength). Detailed protocols and test reports are provided in Section 18 of the submission.

	Test Results				
Test	Proposed Device: The Zimmer Nevel Total Elbow	Predicate Device: The Coonrad/Morrey Total Elbow (K001918, K053189)			
Fatigue Testing: Humeral Stem (Plasma Spray Region)	Five Zimmer Nexel Total Elbow humeral stems achieved (b)(4) load cycles at (b)(without fracture.	The predicate device humeral component was fatigue tested at various load levels, and the fatigue strength at $(b)(4)$ cycles was determined to be (b)			
Fatigue Testing: Ulnar Stem (Mid-Stem Region)	Five Zimmer Nexel Total Elbow ulnar stems achieved (b)(4) load cycles in the mid-stem region at a load of (b)(- without fracture.	The predicate device ulnar component was fatigue tested at various load levels, and the fatigue strength at $(b)(4)$ load cycles in the mid-stem region was determined to be $(b)(4)$			

Fatigue Testing: Ulnar Stem (Plasma Spray Region)	Five Zimmer Nexel Total Elbow ulnar stems achieved $\frac{1}{(b)}$ load cycles in the plasma spray region at a load of (b) without fracture.	The predicate device ulnar component was fatigue tested at various load levels, and the fatigue strength at $(b)(4)$ load cycles in the plasma spray region was determined to be $(b)(4)$
Wear Testing	The Zimmer Nexel Total Elbow mean gravimetric wear rate (b)(4) =: at (b)(4) load cycles) was (b) less than that of the predicate.	The predicate device mean gravimetric wear rate was (b)(4) load cycles.
Durability Testing	The Zimmer Nexel Total Elbow device achieved $(b)(4)$ outs at an equivalent of $(b)(-)$ weight-in-hand.	The run out load for the predicate device was found to be equivalent to (h)(4) weight-in-hand.
Modular Connection Fatigue Testing	The median fatigue strength of the non- articulating, mechanically locked, modular Zimmer Nexel Total Elbow implant components was no less than (b)(4.	Not Applicable: The Zimmer Coonrad/Morrey Total Elbow does not have non-articulating, mechanically-locked, modular implant components comparable to those of the proposed device.


The testing confirms that the Nexel Total Elbow is capable of withstanding expected in-vivo loading and is substantially equivalent to predicate elbow systems.

The sponsor also performed bench testing to support a claim regarding delamination resistance of the Vivacit-E UHMWPE bearing material. Michael Kasser, Ph.D., polymers expert in DOD, acknowledged that this claim was similar to previous claims made for this material (in other joints) and appeared acceptable as provided. Casey Hanley, Ph.D., reviewer in JFOB has reviewed similar Zimmer claims for this material and also confirmed that it appeared acceptable.

The claim reads as follows:

Proposed Claim

Vivacit-E Vitamin E Highly Crosslinked Polyethylene exhibits improved delamination resistance compared to conventional polyethylene under aged conditions (per ASTM F2003). The results of *in vitro* delamination tests have not been shown to correlate with clinical delamination mechanisms.



11



Testing supports this claim. See page 207 in Section 18 and the test report 26_ZRR_WA_2580_12 in Section 18b, for details.

Reviewer Comment: The Nexel Total Elbow has been evaluated for humeral and ulnar stem fatigue strength, wear resistance, durability, and modular connection fatigue strength. The testing demonstrates that the Nexel Total Elbow is capable of withstanding expected in-vivo loading and is substantially equivalent to predicate elbow systems. In addition, testing supports the sponsor's proposed claim regarding delamination resistance. Performance testing is adequate. No additional testing is necessary.

XI. Substantial Equivalence Discussion

The materials, design, sizes, indications, intended use, sterilization, and packaging, for the subject device are similar and/or identical to legally marketed predicate devices. The sponsor provided, via interactive review, additional information including EtO residual testing results, details on colorants used in their instrumentation, updated labeling to include appropriate MR statements, additional discussion regarding the differences in range of motion between their device and the predicate device, edits to all relevant parts of the submission to replace the device descriptive term "semi-constrained" with the appropriate term "constrained", and revisions to the 510(k) Summary.

The sponsor has provided adequate characterization in terms of device design, sizes, material, intended used and testing to demonstrate that the subject device components are equivalent to marketed predicate devices. <u>Therefore, I recommend SE.</u>

		Yes	No	
1.	Same Indication Statement?	X	1	If YES = Go To 3
2.	Do Differences Alter The Effect Or Raise New Issues of Safety Or Effectiveness?	ł		If YES = Stop NSE
3.	Same Technological Characteristics?	X	1	If YES = Go To 5
4.	Could The New Characteristics Affect Safety Or Effectiveness?			If YES = Go To 6
5.	Descriptive Characteristics Precise Enough?		X	If NO = Go To 8 If YES = Stop SE
6.	New Types Of Safety Or Effectiveness Questions?			If YES = Stop NSE
7.	Accepted Scientific Methods Exist?		1	If NO = Stop NSE
8. 1	Performance Data Available?	X	1	If NO = Request Data
9. Data Demonstrate Equivalence?		X	1	Final Decision: <u>SE</u>

Note: See

http://eroom.fda.gov/eRoomReg/Files/CDRH3/CDRHPremarketNotification510kProgram/0_4148/FLOWCHART %20DECISION%20TREE%20.DOC for Flowchart to assist in decision-making process. Please complete the following table and answer the corresponding questions. "Yes" responses to questions 2, 4, 6, and 9, and every "no" response requires an explanation.

Questions 2, 4, and 6 are not applicable.

- 1. Explain how the new indication differs from the predicate device's indication:
- 2. Explain why there is or is not a new effect or safety or effectiveness issue:
- 3. Describe the new technological characteristics:
- Explain how new characteristics could or could not affect safety or effectiveness:
- 5. Explain how descriptive characteristics are not precise enough: Testing, to demonstrate that the design, strength, and compatibility of the subject components are adequate to perform as intended and are comparable to predicate devices, and support proposed wear claims, needs to be provided.
- Explain new types of safety or effectiveness question(s) raised or why the question(s) are not new:
- 7. Explain why existing scientific methods can not be used:
- 8. Explain what performance data is needed: Bench test data (fatigue, wear, modular connections) to demonstrate that the subject device compares favorably to predicate devices. Testing also needs to support proposed claims regarding delamination resistance of the Vivacit-E UHMWPE.
- 9. Explain how the performance data demonstrates that the device is or is not substantially equivalent: Bench test data provided in this submission, including fatigue strength, wear, durability, and modular connection testing demonstrated that the subject device performs as well as or better than predicate devices. Testing also supports proposed claims regarding delamination resistance of the Vivacit-E UHMWPE.

XII. Deficiencies

None.

XIII. Contact History

12/28/12 – Sponsor was notified via email of the RTAX (Refuse to Accept – transition period) status of their submission. A copy of the RTA checklist was provided, that included comments on those items noted to be deficient. See submission file for details.

1/15/13 – Phone conversation with Ms. Joanna Surma discussing the timeline for Zimmer's submittal of EO residual data for FDA review. See submission file (1/16/13 email from Ms. Surma and attachment with RTA response dated 1/17/13) for details.

1/16/13 – Sponsor responded via email to the deficiencies cited in the RTA checklist. All items were adequately addressed. See amended submission file for details.

2/11/13 - Sent an email to the sponsor requesting the EO residual test results.

2/13/13 – Automated email sent to sponsor notifying them that a substantive review of their submission had been completed and that FDA would proceed interactively with them to resolve any remaining deficiencies.

2/14/13 - Sponsor provided EO residual test results via email.

2/22/13 – Phone conversation with the Ms. Surma and Jason Heckaman to discuss outstanding issues to be addressed interactively. General issues were forwarded to the sponsor via email prior to the telecon. They included the following:

- Need for MR statement in the package insert
- Identify any color additives that may be included in any device instrumentation.
- Update 510(k Summary to include correct regulation/classification name, procode name. Semiconstrained versus constrained – revise.
 - Classification: 21 CFR 888.3150 Elbow joint metal/polymer constrained cemented prosthesis.
 - o Product Code: JDC Prosthesis, Elbow, Constrained, Cemented
- Nitriding location: on interior surface of ulnar eye or just outer surface, or both?
- Variations in ROM from predicate; discussion of differences and the effect, if any.

2/27/13 - Sponsor responded via email to the requests made during the 2/22/13 teleconference. They provided all the requested information. See appropriate sections of this review memo for details about the responses.

2/28/13 – Sponsor provided an updated 510(k) Summary via email. See Reviewer Comment in the Labeling section of this review memo for details.

All email correspondence and attachments have been added to the submission file. See submission file, as amended, for details.

XIV.Recommendation

Regulation Number: 21 CFR 888.3150 Regulation Name: Elbow joint metal/polymer constrained cemented prosthesis Regulatory Class: Class II Product Code(s): JDC

I recommend this submission, as revised by emails submitted during interactive review, be found SE. An SE Letter should be sent to the sponsor.

Digital Signature Concurrence Table			
Reviewer Sign-Off	Peter G. Allen 2013.03.08 16:09:20 -05'00' DA		
Branch Chief Sign-Off			
Division Sign-Off			

14

From:	loanna I. Surma < lo	anna.Surma@zimmer	com>	
Sent:	Thursday, Eebruary	Super or modeling from a la de trais or sinds descada		
To:	Allen, Peter			
Subject:		will Proceed via Intera	ctive Review	
Attachments:		10(k) Summary.pdf 👉		
Importance:	High			
Hi Pete,				
I noticed that I provided th	he wrong Attachment 2. Plea	ase replace the one pr	eviously sent with t	this one.
Thanks,				
Joanna				
Joanna L. Surma - MS, RA	c			
The set of helpersphere is the first state that the set of the first state of the set of				
Associate Project Manage Regulatory Affairs	I			
Zimmer, Inc.				
(o) 574.371.1642				
(c) 574.527.2361				
From: Joanna L. Surma		n, (/) (/ ================================		
Sent: Wednesday, Februa Tot, 'Allen, Peter'	iry 27, 2013 1:48 PM			
	will Proceed via Interactive	Review		
Importance: High				
Hi Pete,				
Thanks for the call last we	ek. Below, are the clarificati	ons requested.		
 MR Statements in insert. 	the Package Insert: See Atta	achment 1 (page 4) fo	r a revision of the N	lexel Total Elbow package
	ditives in non-metallic instru with the same colorant add		n part numbers and	K-numbers of legally
Non-Metallic Part Number	Description	Non-Metallic Material	Color Additive?	Part Number and K-Numbe of Legally Marketed Device with the Same Color Additive
00-8401-028-00	Ulnar Stem Inserter	Radel	none	NA
00-8491-088-00	Humeral Screw Installation Driver Bit	LDPE	none	NA



P.O. Box 708 Warsaw, IN 46581-0708 574 267-6131

ì

510(k) Summary of Safety and Effectiveness

Sponsor:	Zimmer, Inc. P.O. Box 708 Warsaw, IN 46581-0708				
Contact Person:	Joanna L. Surma Associate Project Manager, Regulatory Affairs Telephone: (574) 371-1642 Fax: (574) 372-4605				
Date:	13 December 2012				
Trade Name:	Zimmer Nexel Total Elbow				
Common Name:	Total Elbow Prosthesis				
Classification	21 CFR § 888.3150 Elbow joint metal/polymer constrained cemented prosthesis				
Product Code:	JDC - Prosthesis, Elbow, Constrained, Cemented				
Predicate Device:	Coonrad/Morrey Total Elbow, manufactured by Zimmer, <u>K001989</u> , cleared 25 July 2000				
	Coonrad/Morrey Total Elbow, manufactured by Zimmer, <u>K053189</u> , cleared 9 December 2005				
Device Description:	The Zimmer Nexel total Elbow is a total elbow prosthesis designed for use with bone cement. It is available in multiple sizes and in right and left configurations.				
	How the Device Functions: The Zimmer Nexel Total Elbow is an implant designed to replace the articulating surfaces of and restore motion to the humeroulnar (elbow) joint. The implant is a constrained device assembly and consists of the following components: humeral component, ulnar component, humeral bearing-A, ulnar bearings-B, axle pin and humeral screws.				
х.	The Following Scientific Concepts, Design Features and Physical Properties form the Basis for the Zimmer Nexel Total Elbow: The humeral component has a humeral stem designed to be implanted with bone cement into the patient's humeral medullary canal, an anterior flange designed to accept a				

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

×.

510(k) Summary of Safety and Effectiveness – Zimmer Nexel Total Elbow

bone graft and limit torsional and posterior migration, a humeral yoke with rounded corners to avoid the creation of stress risers within the medial and lateral humeral supracondylar columns, and plasma spray region to enhance fixation to bone cement within the medullary canal, and to improve fatigue strength. The ulnar component has an ulnar stem designed to be implanted with bone cement into the patient's ulnar medullary canal, an ulnar eye that is both highly polished and nitrogen-enriched to limit wear of the apposing polymer bearings, and plasma spray region to enhance fixation to bone cement within the medullary canal. Bearings A and B are designed to broadly distribute joint reaction forces.

Materials Used: The humeral and ulnar components are made of a titanium alloy, the bearings A and B are made of Vitamin E highly cross-linked polyethylene, and the humeral screws and axle pin are made of a cobalt-chromium-molybdenum alloy.

Comparison to the Predicate: The proposed device (Zimmer Nexel Total Elbow) and the predicate device (Coonrad/Morrey Total Elbow) have the same intended use and similar indications for use. The proposed device humeral and ulnar components are very similar in terms of materials used and design/dimensions to the predicate device humeral and ulnar components. The primary differences between the proposed and predicate devices are that the proposed device bearing components are made of Vitamin E highly cross-linked polyethylene (the predicate device bearing components are made of ultra-high molecular-weight polyethylene), and the proposed device has bearings that articulate on both the outer and inner diameters of the ulnar eye (the predicate device has bearings that articulate on the inner diameter of the ulnar eye). Additionally, the proposed device uses screws and a 1-part axle pin to fix the bearings in place, while the predicate device uses a 2-part (snap-fit) axle pin.

Intended Use:

Indications for use include:

- Elbow joint destruction which significantly compromises the activities of daily living
- Post-traumatic lesions or bone loss contributing to elbow instability
- Ankylosed joints, especially in cases of bilateral ankylosis from causes other than active sepsis
- Advanced rheumatoid, post-traumatic, or degenerative arthritis with incapacitating pain
- Instability or loss of motion when the degree of joint or soft tissue damage precludes reliable osteosynthesis
- Acute comminuted articular fracture of the elbow joint surfaces that precludes less radical procedures, including 13-C3 fractures of the distal humerus
- Revision arthroplasty

Caution: This device is intended for cemented use only.

510(k) Summary of Safety and Effectiveness - Zimmer Nexel Total Elbow

Comparison to Predicate Device:	The Zimmer Nexel Total Elbow is substantially equivalent to the predicate device in terms of form and function. The Zimmer Nexel Total Elbow and the predicate device share similar intended uses and indications for use.
Performance Data:	 Non-Clinical Performance Testing Conducted: Stem Fatigue Testing Wear Testing Durability Testing Modular Connection Fatigue Testing
	Non-Clinical Performance Testing Conclusions:

Non-clinical testing demonstrated that the New Zimmer Total Elbow meets performance requirements as defined by Design Control activities and is substantially equivalent to the predicate device in terms of safety and efficacy.

In this case, clinical data and conclusions were not needed to demonstrate substantial equivalence.

lFrom:	Joanna L. Surma <joanna.surma@zimmer.com></joanna.surma@zimmer.com>			
Sent:	Wednesday, February 27, 2013 1:48 PM			
То:	Allen, Peter			
Subject:	K123862 Review will Proceed via Interactive Review			
Attachments: 🥧	Attachment 113c. Package Insert.pdf;/Attachment.3;- 18b.7 ZRR_WA_2476_11.pdf; Attachment 2:-05. 510(k) Summary.pdf; Attachment 5Device Descriptionpdf; Attachment 6 ** RTA Response - 02 - Zimmer Responses.pdf;-Attachment 7 - RTA Response - 12a. Substantial Equivalence Discussion.pdf;*Attachment 4*- Executive Summary.pdf			
Importance:	High			

Allen, Peter

Hi Pete,

Thanks for the call last week. Below, are the clarifications requested.

- 1. MR Statements in the Package Insert: See Attachment 1 (page 4) for a revision of the Nexel Total Elbow package insert.
- 2. List of colorant additives in non-metallic instrumentation along with part numbers and K-numbers of legally marketed devices with the same colorant additives:

Non-Metallic Part Number	Description	Non-Metallic Material	Color Additive?	Part Number and K-Number of Legally Marketed Device with the Same Color Additive
00-8401-028-00	Ulnar Stem Inserter	Radel	none	NA
00-8401-088-00	Humeral Screw Installation Driver Bit	LDPE	none	NA
00-8401-089-00	Humeral Screw Removal Driver	Silicone	Yes	00-2360-087-00 K040593
00-8401-080-00	Torque-Limiting Humeral Screw Driver	Ultem	Yes	00-5971-050-09 K933785

- 3. Update 510(k) summary to include correct regulation/classification name and product code name: See Attachment 2 for an updated 510(k) summary.
- 4. Nitriding location: On all surfaces of ulnar eye. Discussion point closed during phone call.
- 5. Discussion of ROM:

The Nexel Total Elbow's allowable flexion arc from (b)(4) and varus/valgus laxity of (b)(4) comply well with the 0-146° flexion arc and +/-3°-5° range of varus/valgus angulation of a normal, healthy elbow reported in the referenced studies summarized in Zimmer's BIOMECHANICS RATIONALE FOR THE LOADING OF TOTAL ELBOW PROSTHESIS DURING IN VITRO SIMULATIONS document (included in the initial K123862 submission and also provided here as Attachment 3). While the predicate C/M implant assembly allows (b)(4) (depending on the

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

size configuration) flexion and up to (b)(of v/v laxity, the normal elbow doesn't allow (due to soft and hard tissue limitations) this much flexion and functions normally within the (b)(4 - v/v angulation range provided by both the Nexel and the C/M.

The Nexel Total Elbow was designed to provide at least as much of the necessary flexion and v/v laxity as required of a normal, healthy elbow while maximizing the thickness and optimizing the location and form of the Vitamin-E UHMWPE Bearings to provide improved wear resistance and longer-term load-bearing capability as compared to the predicate device. Thus, any flexion and v/v angulation deemed unnecessary for a full healthy elbow range of motion was traded for optimization of device performance; this is further substantiated by the results of the head-to-head Wear and Durability tests reported in the initial K123862 submission, wherein both (Nexel and predicate) devices were subjected to high-loading over large flexion arcs with induced v/v angulation to simulate clinically-relevant, yet aggressive, long-term use in-vivo.

- 6. Although, we believe the Nexel Total Elbow is a semi-constrained device in the engineering sense (that it is not a rigid hinge-like device fully-constrained to only rotation about a single axis; that it, in fact, allows for an additional range of varus/valgus motion about that hinge), we agree that this engineering understanding of "semi-constrained" still falls under the FDA definition of "constrained" and have changed the following submission pages accordingly:
 - Initial K123862 submission page 16 see page 1 of Attachment 2
 - Initial K123862 submission page 56 see page 1 of Attachment 4
 - Initial K123862 submission page 61 see page 1 of Attachment 5
 - Initial K123862 submission page 74 see page 14 of Attachment 5
 - RTA Response page 4 see page 2 of Attachment 6
 - RTA Response page 27 see page 9 of Attachment 7
 - Initial K123862 submission page 57 see page 2 of Attachment 4

Please let me know if you require any further clarifications.

Sincerely, Joanna

Joanna L. Surma - MS, RAC Associate Project Manager Regulatory Affairs Zimmer, Inc. (o) 574.371.1642 (c) 574.527.2361

From: Allen, Peter [<u>mailto:Peter.Allen@fda.hhs.gov</u>] Sent: Friday, February 22, 2013 1:42 PM To: Joanna L. Surma Subject: RE: K123862 review will proceed via interactive review

Joanna,

Very generally, here are the issues:

- Need for MR statement in the package insert
- Colorant additives that may be included in any instrumentation
- Update 510(k Summary to include correct regulation/classification name, procode name. Semi-constrained versus constrained revise.

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017

- c Classification: 21 CFR 888.3150 Elbow joint metal/polymer constrained cemented prosthesis.
- o Product Code: JDC Prosthesis, Elbow, Constrained, Cemented
- Nitriding location: on interior surface of ulnar eye or just outer surface, or both?
- Variations in ROM from predicate; discussion of differences and the effect, if any.

Pete

From: Joanna L. Surma [mailto:Joanna.Surma@zimmer.com] Sent: Friday, February 22, 2013 10:28 AM To: Allen, Peter Subject: RE: K123862 review will proceed via interactive review

3pm EST should work for me. Could you please let me know the topics of discussion, so I can make sure to have the appropriate subject-matter experts in the conference?

Thanks, Joanna

Joanna L. Surma - MS, RAC Associate Project Manager Regulatory Affairs Zimmer, Inc. (o) 574.371.1642 (c) 574.527.2361

From: Allen, Peter [mailto:Peter.Allen@fda.hhs.gov]
Sent: Friday, February 22, 2013 10:25 AM
To: Joanna L. Surma
Subject: RE: K123862 review will proceed via interactive review

Joanna,

I have a couple of items I'd like for you to address for me before I SE this submission. I will attempt to call you later this afternoon to discuss (~3pm) if my schedule permits (I'm in meetings most of the day today), otherwise I will send an email to you and we can discuss next week.

Pete Allen 301-796-6402

From: Joanna L. Surma [mailto:Joanna.Surma@zimmer.com] Sent: Monday, February 18, 2013 6:09 PM To: Allen, Peter Subject: RE: K123862 review will proceed via interactive review

Hello Mr. Allen,

Please let me know how you'd like to proceed with the interactive review.

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

Thanks, and I look forward to hearing from you soon,

Joanna

Joanna L. Surma - MS, RAC Associate Project Manager Regulatory Affairs Zimmer, Inc. (o) 574.371.1642 (c) 574.527.2361

From: Peter Allen [mailto:peter.allen@fda.hhs.gov] Sent: Wednesday, February 13, 2013-3:52 PM To: Joanna L. Surma; Jason Heckaman Cc: Peter Allen Subject: K123862 review will proceed via interactive review

February 13, 2013 Substantive Review Notification - Proceed Interactively

FDA has completed a substantive review of your submission K123862 and has determined that we will work with you to resolve outstanding deficiencies through an Interactive Review process. Any outstanding deficiencies will be provided to you during the Interactive Review process.

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017 Attachment / Bar code 87-5204-077-99 Rearly CE CE CE

Cé ECREP Zimmer, Inc Zimmer U.K. Ltd. 1800 West Center Street 9 Lancaster Place Warsaw, Indiana 46580 South Marston Park zimmer USA Swindon, SN3 4FP, 1-800-348-2759 (US only) +1-574-372-4999 UK

English

The CE mark is valid only if it is also printed on the product label

ZIMMER® NEXEL ® TOTAL ELBOW

Before using a product placed on the market by Zimmer, the operating Surgeon should study carefully the following recommondations, warnings and instructions, as well as the available product-specific information (e.g. product threature, written surgical technique). Zimmer is not liable for complications arising from the use of this device outside of its indicated uses, not adgined with the surgical technique, product selection or similar matters outside the control of Zimmer.

DESCRIPTION

DESCRIPTION This denote is a total above prostheses designed for use with bone cement. It is available in sizes 4, 5 and 8, in left and right configurations. The uner and numeral implant components are manufactured from *Twanium*® T+6A-4V alloy. The uner implant component has a produe boaring of T+6A-4V plasma spray and is curred to faoilates implantation. The humeral implant component has a produe boaring of T+6A-4V plasma spray and has an anterior flange to accommodate a bone graft. The inhage prima and servers are manufactured from *Zimatoy*® Co-Cr-Mo alloy. Mamin Einghly cross-Inked ultra-high molecular weight polyathylene (*Vwoch E*®) bearings prevent metal-to-metal antioutating contact.

INDICATIONS

- Indications for use include

- Instants for use proutee Ellow pain destruction which significantly comprises the activities of daily living Post-traumatic lessons or bone loak contributing to elbow instability Ankytosed joints, especially in cases of bilatoria inkytosis from causes other han solve sepses Advanced theumatoid, post-tauramitic, or degenerative artification apacitating pan instability or loss of motion when the degree of joint or soft insue damage predudes reliable
- ostaceynthesis
- Acute comminuted articular fracture of the elbow joint surfaces that produdes less radical
 procedures, including 13-C3 fractures of the distal humans
 Revision arthroptacty

CALITION: This device is intended for camenied use only

- CONTRAINDICATIONS Use of the Zimmer Nexe/Total Elbow is contraindicated in patients with
- Currently active, or history of repeated, local infection at the surgical arts
 Paralysis or dysfunctional neuropathy involving the elbow joint

- Significant positive hand dystunction Excessive scamp of the ston or soft basis that could prevent adequate soft tissue coverage Daily acrivities that would subject the device to significant stress (i.e., heavy labor, torsional stress, and/or competitive sports)

- Relative contrandications include Distant foor of infection (e.g. genitourinary, pulmonary, skin [chronic leaces or uccatations] or other *stasy*, *1*2 cases of detent infection, the loss of intection should be treated prior to, during and effer-ments or the standard statements of the standard statement of the standard statement of the standard statement of the statement
- surgery. Ancient prior seosis

WARNINGS

- WARNINGS
 This device is for angle patient use only. Do not reuse. Rouse of a single use device that has come in contact with blood, bone, tassue or other body flucts may lead to patient or user injury. Possible miss associated with reuse of a engle use device include, but are not limited to, mechanical failure and transmission of infectious egents.
 To properly implant this device the surgical technique should be consulted and carefully followed.
 Carmenting technique is extremely important. The medulary canal should be copyoutly implant the device in the interval limit device of the surgical technique should be consulted and carefully followed.
 Carmenting technique is extremely important. The medulary canal should be copyoutly implant to remove blood, fat, and bern deating and then throughly dired. A cement derivery system should be used to carrient the humanal and ulnar implant components.
 A sing mechanical Right for the implant within the humanal and ulnar medulary canels may help minimize isocareing of the implant upplant within the transmission.

- A anug mechanical fit of the implant within the humanal and ultrar meduliary canels may help minin iconsening of the implant components. The smount of bone removed from the ultra should be sufficient to permit full elbow motion on the operating table when the appropriate provisional prositions is fully inserted. The surgeon should check for full range of elbow motion at appropriate times during surgery.

PRECAUTIONS Proper handli

- PRECAUTIONS
 Proper handling of this implant is important. Contouring (bending) of the humanal or uthan elems should be avoided. Any attention may produce defects that could become focal points for excessive affects leading to implant failure. To prevent talse infection of lotal point mplants, many surgeons advice the use of antibiotic prophytaxs before and after dental manipulation, endoscopic examinations and other minor surgical procedures. Patients with eignificant lower activenity disability who require walleng and a me leas amenable to total eabore attrophasy tion patients with advanced and precommany upper extremity involvement. The *Zimmer Nexet* Total Elbow has not been evaluated for safety and compatibility in the MR economical.

- environment
 The Zammer Name/ Total Etocw has not been evaluated for heating or migration in the WR
- ADVERSE EFFECTS

- ADVERSE EFFECTS
 The following advance effects have been reported
 Early or late component locating
 Implant failure or fracture
 Infection
 Nerve injury
 Thoops and ation
 Cortical performances and fractures
 Cortical performances and fractures
 Advance lister engones and inflammatory reaction to the implant material or wear debns, including
 met at semicilistic
- metal sensitivity Osteolysis and/or loosening due to wear debre .

STERILITY

The Jonner Alexel Total Elbow mptants are provided stanie (the humeral and unar implant components, and the humeral acrews are stanized by gamma insiderion, indicated by the <u>BIERIE R</u> symbol on the isobeling; and, the activuition is 1 – installing the bearings and the auto pin - is stanized by environments and the standard by the <u>BIERIELE FO</u> symbol on the labeling) and romain stante until the expration date displayed on the package, as long as the package integrity has not been wolkled. Inspect each package prior to use and do not use the implant component if any seal or cavity is damaged or breached or if the expiration date has been expeeded. Once opened, the implant component must be used, discarded or re-steritized

STERILIZATION INSTRUCTIONS

These sterilization instructions are consistent with ANSI/AMI//SO standards and guidelines. They should be used for sterile items that were opened but unuesd • Do not reuse implants

- Do not reuse implants In the event of readvertant loss of startility while propering for aurgery, all metal implant components may be re-startized only once for immediate use. The is subject to the exceptions iside below: <u>DO NOT RE-STERILZE</u> Implant components that have been contaminated with body flucts or biologic debris, or that have been previously implanted Implant components ontaining *Visiot-E* polyethylene Any plastic, loarn or polyethylene Do not use the original plastic cantes or lists for re-startization. Single implant components should be poolspaced in a metical gradue destination pound in the commended . packaged in a medical grade stantization pouch or wrap that conforms to the recommended specifications for staam stantization provided in the table under *Recommended Startization/Re*
- sterilization Specifizations below. Ensure that the pouch or wrap as large enough to contain the implant component without stressing its seals or tearing Implant components must be disassembled and sterilized separately to minimize potential bio-burden.
- Impart componentiate or trace and trace and trace and trace and trace and trace and trace or trace or

.

Recommended Sterilization/Re-sterilization Specifications Follow the sterilizer menufacturer's instructions for loading patterns and selection of elerilization parameters. Drying times vary according to load size and should be increased for larger loads

For Single, All Metal Implante

Sieam Stenizman				
Тур∎	Temperature	Exposure Time	Minimum Dry Time	
Gravity Displacement	121"C (250"F)	30 minutes	15 minules	
Gravity Displacement	132*C (270*F)	15 minules		
UK Pre-vacuum/Pulsating Vacuum ¹	134°C (273°F)	3 minules		
Pre-vacuum/Puisaling Vacuum	132°C (270°F)	4 minutes		

1 This cycle is not for use in the United States

Please contact Zimmer at the following number if you have additional questions. In the USA, call 1-800-348-2759 For calls outside the USA, call the local international access code +1-574-267-5131.

- STERILIZATION PRECAUTIONS
 Appressive descring with detergents and brushes may damage special features of the implants, such as porcus costings. Also, certain detergents may be difficult to mas off polymer items, especially those made for iterum and thermum alloys can form could layers from steam boiler beatment chemicats or detergent readures. While these oxides are brocompatible, they can obligate iterate identification effortings and elempings.

PATIENT COUNSELING INFORMATION Complications analysis in proclimatic implants are more likely to occur in patients, with unrealistic functional expectations, heavy patients, physically active patients, and/or with patients who fail to follow through with the required netabilitation program. Physical soluty can result in locesning, wear, and/or fracture of the implant. The proceedness implant patient must be counseled about the capabilities of the implant and the implact the proceedness implant patient intuits be counseled about the capabilities of the implant and the implant in the proceedness implant patient to couparional and sports activates, and about the possibility that the implant in its components may wear out. If all or need to be replaced the implant may regilitable the rest of the patients if to or any particular length of time. Because prosthetic implants are not as group, reliable, or durable as natural, healthy biscueschores, all such devices may need to be replaced at group point. at some point

Note The proceduress will not reactors function to the level expected with a normal healthy joint, and has lunctional limitations. Exclessive muscular activity such as pounding and or carrying excessive loads must be exclosed or reducted. The patient must not thit mays than one pound (~0.5 kg) during the first inner postopersities months; and thereafter, not more than five pounds (~2.5 kg) with the operated arm in addition, the rests and complications listed under ADVERSE EFFECTS must be explained to and discussed with the patient

CAUTION: This device is intended for comented use only.

.

· .



Date:

Trade Name:

Classification

Product Code:

Predicate Device:

Device Description:

Common Name:

P.O. Box 708 Warsaw, IN 46581-0708 574 267-6131

510(k) Summary of Safety and Effectiveness

Zimmer, Inc.		
P.O. Box 708		
Warsaw, IN 46581-0708		

anna L. Surma		
Associate Project Manager, Regulatory Affairs		
lephone: (574) 371-1642		
x: (574) 372-4605		
1		

13 December 2012

Zimmer Nexel Total Elbow

Total Elbow Prosthesis

21 CFR § 888.3150 Elbow joint metal/polymer constrained cemented prosthesis

JDC – Prosthesis, Elbow, Constrained, Cemented

Coonrad/Morrey Total Elbow, manufactured by Zimmer, K001989, cleared 25 July 2000

Coonrad/Morrey Total Elbow, manufactured by Zimmer, K053189, cleared 9 December 2005

The Zimmer Nexel total Elbow is a total elbow prosthesis designed for use with bone cement. It is available in multiple sizes and in right and left configurations.

How the Device Functions: The Zimmer Nexel Total Elbow is an implant designed to replace the articulating surfaces of and restore motion to the humeroulnar (elbow) joint. The implant is a semi-constrained device assembly and consists of the following components: humeral component, ulnar component, humeral bearing-A, ulnar bearings-B, axle pin and humeral screws.

The Following Scientific Concepts, Design Features and Physical Properties form the Basis for the Zimmer Nexel Total Elbow: The humeral component has a humeral stem designed to be implanted with bone cement into the patient's humeral medullary canal, an anterior flange designed to accept a

510(k) Summary of Safety and Effectiveness - Zimmer Nexel Total Elbow

bone graft and limit torsional and posterior migration, a humeral yoke with rounded corners to avoid the creation of stress risers within the medial and lateral humeral supracondylar columns, and plasma spray region to enhance fixation to bone cement within the medullary canal, and to improve fatigue strength. The ulnar component has an ulnar stem designed to be implanted with bone cement into the patient's ulnar medullary canal, an ulnar eye that is both highly polished and nitrogen-enriched to limit wear of the apposing polymer bearings, and plasma spray region to enhance fixation to bone cement within the medullary canal. Bearings A and B are designed to broadly distribute joint reaction forces.

Materials Used: The humeral and ulnar components are made of a titanium alloy, the bearings A and B are made of Vitamin E highly cross-linked polyethylene, and the humeral screws and axle pin are made of a cobalt-chromium-molybdenum alloy.

Comparison to the Predicate: The proposed device (Zimmer Nexel Total Elbow) and the predicate device (Coonrad/Morrey Total Elbow) have the same intended use and similar indications for use. The proposed device humeral and ulnar components are very similar in terms of materials used and design/dimensions to the predicate device humeral and ulnar components. The primary differences between the proposed and predicate devices are that the proposed device bearing components are made of Vitamin E highly cross-linked polyethylene (the predicate device bearing components are made of ultra-high molecular-weight polyethylene), and the proposed device has bearings that articulate on both the outer and inner diameters of the ulnar eye (the predicate device has bearings that articulate on the inner diameter of the ulnar eye). Additionally, the proposed device uses screws and a 1-part axle pin to fix the bearings in place, while the predicate device uses a 2-part (snap-fit) axle pin.

Intended Use:

Indications for use include:

- Elbow joint destruction which significantly compromises the activities of daily living
- Post-traumatic lesions or bone loss contributing to elbow instability
- Ankylosed joints, especially in cases of bilateral ankylosis from causes other than active sepsis
- Advanced rheumatoid, post-traumatic, or degenerative arthritis with incapacitating pain
- Instability or loss of motion when the degree of joint or soft tissue damage precludes reliable osteosynthesis
- Acute comminuted articular fracture of the elbow joint surfaces that precludes less radical
 procedures, including 13-C3 fractures of the distal humerus
- Revision arthroplasty

Caution: This device is intended for cemented use only.

510(k) Summary of Safety and Effectiveness - Zimmer Nexel Total Elbow

Comparison to Predicate Device: The Zimmer Nexel Total Elbow is substantially equivalent to the predicate device in terms of form and function. The Zimmer Nexel Total Elbow and the predicate device share similar intended uses and indications for use.

Performance Data:

Non-Clinical Performance Testing Conducted:

- Stem Fatigue Testing
- Wear Testing
- Durability Testing
- Modular Connection Fatigue Testing

Non-Clinical Performance Testing Conclusions:

Non-clinical testing demonstrated that the New Zimmer Total Elbow meets performance requirements as defined by Design Control activities and is substantially equivalent to the predicate device in terms of safety and efficacy.


In this case, clinical data and conclusions were not needed to demonstrate substantial equivalence.

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017 2/27/13 e mái



Page 1 of 29





.

2/27/13 email



Executive Summary Premarket Notification

Device Name

Zimmer Nexel Total Elbow


Reason for the Submission


Zimmer is seeking clearance to market the Zimmer Nexel Total Elbow, which has not previously been marketed.

Device Description

Overview

The Zimmer Nexel Total Elbow is designed to replace the articulating surfaces of and restore motion to the humeroulnar (elbow) joint. The implant is a constrained device assembly and consists of a *humeral component*, an *ulnar component*, a *humeral bearing-A*, 2 *ulnar bearings-B*, an *axle pin* and 2 *humeral screws*.





A detailed description of the proposed device is presented in <u>Section 11</u>, <u>Device</u> <u>Description</u>.



Indications for Use

Indications for use include:

- Elbow joint destruction which significantly compromises the activities of daily living
- Post-traumatic lesions or bone loss contributing to elbow instability
- Ankylosed joints, especially in cases of bilateral ankylosis from causes other than active sepsis
- Advanced rheumatoid, post-traumatic, or degenerative arthritis with incapacitating pain
- Instability or loss of motion when the degree of joint or soft tissue damage precludes reliable osteosynthesis
- Acute comminuted articular fracture of the elbow joint surfaces that precludes less radical procedures, including 13-C3 fractures of the distal humerus
- Revision arthroplasty

CAUTION: This device is intended for cemented use only.

Submission Overview

The table below provides a comparison between the proposed and predicate devices. Additional information is presented in <u>Section 12</u>, <u>Substantial Equivalence</u> <u>Discussion</u>.

		Proposed Device: The Zimmer Nexel Total Elbow	Predicate Device: The Coonrad/Morrey Total Elbow (K001918, K053189)
	Stem with Triangular Cross-Section	Х	X
글로	Straight Stem	X	<u>x</u>
Humeral Component	Anterior Flange	x	x
	Titanium Alloy	Х	X
± 3	Titanium Beads on Distal Region		x
	Titanium Plasma Spray on Distal Region	X	
	Stem with Quadrangular Cross-Section	<u>X</u>	x
LT L	Curved Stem	X	x
Ulnar Component	Nitrogen Enriched Surface Hardening of Ulnar Eye	х	
Ü	Titanium Alloy	X	x
	Titanium Plasma Spray on Proximal Region	x	x
- e.	2-Part, Hollow Snap-Pin		x
Linkage	Solid Axle Pin	х	
Ľ.	Screws	Х	
	Constrained	x	X
	Terminal Sterilization	х	x
	Titanium-on-Polyethylene Articulation	Х	X



Performance Testing

The following performance testing was conducted on the Zimmer Nexel Total Elbow. A detailed summary and copy of each test report are included in <u>Section 18, Performance</u> <u>Testing – Bench</u>.

Stem Fatigue Testing

The fatigue strength of the Zimmer Nexel Total Elbow humeral component distal stem at (b)(4) load cycles is (b)(-, which is (b)(- stronger than that of the legally marketed Zimmer Coonrad/Morrey Total Elbow humeral component (b)(4) at a comparable point on the distal stem.

The fatigue strength of the Zimmer Nexel Total Elbow ulnar component stem in the mid-stem region at (b)(4) load cycles is at least (b)((the fatigue strength of the legally marketed Zimmer Coonrad/Morrey Total Elbow ulnar component at a comparable point on the stem).

The fatigue strength of the Zimmer Nexel Total Elbow ulnar component stem in the plasma spray (proximal stem) region at (b)(4) load cycles is at least (b)/. I (the fatigue strength of the legally marketed Zimmer Coonrad/Morrey Total Elbow ulnar component at a comparable point on the stem).

These studies demonstrated that the Zimmer Nexel Total Elbow is at least as safe and effective as the predicate device (Zimmer Coonrad/Morrey Total Elbow) in terms of humeral and ulnar stem fatigue strength.

Wear Testing

The Zimmer Nexel Total Elbow mean gravimetric wear rate (b)(4) at (b)(4) load cycles) was (b), less than that of the predicate.

This study demonstrated that the Zimmer Nexel Total Elbow met or exceeded all acceptance criteria and is at least as safe and effective as the predicate device in terms of wear performance.

Durability Testing

The Zimmer Nexel Total Elbow device achieved (b)(4) prunouts at the equivalent of (b)(4) weight-in-hand (b)(4) peak JRF and (b)(4) varus-valgus moment), a load level (b) higher than the load level at which the predicate device ran out.

This study demonstrated that the Zimmer Nexel Total Elbow is at least as safe and effective as the predicate device in terms of durability.



Modular Connection Fatigue Testing

The median fatigue strength of the non-articulating, mechanically locked, modular Zimmer Nexel Total Elbow implant components was no less than (b)(4

This study demonstrated that the Zimmer Nexel Total Elbow is safe and effective in terms of non-articulating, mechanically locked, modular connection fatigue.

Support Research

The following internal (Zimmer) research and testing informed and/or support the aforementioned testing and can be found in their entirety at the end of <u>Section 18</u>, <u>Performance Testing – Bench</u>.

- <u>Zimmer Research Memo ZRM WA 0179 09, Analysis of Retrieved</u> <u>Coonrad/Morrey Elbow Total Arthroplasty Components</u>
- <u>Zimmer Research Memo ZRM_WA_0193_09</u>, Development of the <u>Coonrad/Morrey Elbow Durability Test</u>
- <u>Zimmer Research Report ZRR_WA_1988_09</u>, Literature Review Summary of <u>Elbow Joint Mechanics and Total Elbow Arthroplasty Outcomes</u>
- Zimmer Research Report ZRR_WA_2252_10, Evaluation of the AS/BF Taper Adaptor Component Using Electrochemical Test Methods
- Zimmer Technical Memo 1222.00, Fatigue Strength of Plasma Sprayed Ti-6AL-4V Alloy Material
- <u>Zimmer Research Memo ZRM_WA_0292_11, Summary of Clinical</u> <u>Complications of the Coonrad/Morrey Elbow System</u>
- Zimmer Research Report ZRR WA 2338 11, Vitamin E Highly Crosslinked Polyethylene Shelf Life Study
- Zimmer Research Report ZRR_WI_2441_11, Finite Element Analysis of the Coonrad-Morrey Total Elbow System
- <u>Zimmer Research Report ZRR_WA_2476_11</u>, Biomechanics Rationale for the Loading of Total Elbow Prosthesis During In Vitro Simulations
- <u>Zimmer Research Report ZRR_WA_2481_11, Determination of the</u> <u>Functional Torque Range of the New Zimmer Total Elbow Humeral Screw</u> and <u>Hex Driver</u>
- <u>Zimmer Technical Memo 1217.00, Coonrad/Morrey Elbow Ulnar Stem</u> <u>Fatigue Testing</u>
- <u>Zimmer Technical Memo 1218.00, Coonrad/Morrey Humeral Stem Fatigue</u> <u>Testing</u>
- <u>Memo to DHF File Z07-014, Rationale for Determination of Torque</u> <u>Specification and Torque Range Requirement for the Mechanical Fastening</u> <u>Used in Zimmer New Total Elbow Prosthesis</u>



- <u>Zimmer Research Report ZRR_WA_2646_12, Pin-on-Flat Wear Comparison</u> of the New Zimmer Total Elbow and the Coonrad-Morrey Articulation <u>Couples</u>
- Zimmer Research Report ZRR_WA 2648 12, Finite Element Analysis of Nexel Elbow Ulnar Implant Neck Region
- Zimmer Research Report ZRR_WI_2441_11, Finite Element Analysis of the Coonrad-Morrey Total Elbow System

Records processed under FOIArequest #2016-4653; Released by CDRH on 07/16/2017



2/27/13 email

Traditional 510(k) Premarket Notification

Device Description

Overview

The Zimmer Nexel Total Elbow is an implant designed to replace the articulating surfaces of and restore motion to the humeroulnar (elbow) joint. The implant is a constrained device assembly and consists of the following components, which are described in more detail below:

- Humeral Component
- Ulnar Component
- Humeral Bearing-A
- 2 Ulnar Bearings-B
- <u>Axle Pin</u>
- 2 Humeral Screws








Image 11.2: The Zimmer Nexel Total Elbow, assembled





Image 11.3: The Zimmer Nexel Total Elbow, assembled with Humeral Component translucent



Humeral Component

The *humeral component* is made from wrought *Tivanium* alloy (Ti-6Al-4V) and has the following major design features: a stem, an anterior flange, a yoke and plasma spray regions.



Humeral Component, major design features

The *humeral stem* is intended to be implanted with bone cement into the patient's humeral medullary canal. It has a triangular cross-section, which is designed to prevent rotation of the *humeral component* within the patient's medullary canal; and, it is straight along its proximal-distal axis to conform to the straight shaft of the distal humeral medullary canal.

The *humeral anterior flange* is designed to accept a bone graft, which is wedged between the *anterior flange* and the exterior cortical bone of the humeral shaft. This *anterior flange* and bone graft combination are intended to improve both antirotational stability and posterior migration of the *humeral component*, thereby resisting stem loosening.

The *humeral yoke* has rounded corners (as opposed to angular) on its proximal surface in order to prevent stress risers (that are often associated with bone preparation cuts required for angled corners abutting bone). This rounded *yoke* design is intended to preserve the strength of (and avoid fracture of) the medial and lateral supracondylar columns. The *yoke* walls also house partially threaded holes with Spiralock® technology (tap manufactured by Stanley Black & Decker, Madison Heights, MI), which engage with the *humeral screws*.



The *humeral component* has several regions covered in plasma spray, which is contained within the Plasma Sprayed Titanium-6AI-4V Alloy Coating Material Master File (FDA Ref: MAF-1909, Letter of Authorization can be found in Section 21). The *plasma spray regions* are designed to increase the coefficient of friction between the implant and the bone cement, thereby improving anti-rotational stability and aiding in resistance to movement of the implant with respect to the humerus.

In addition to the aforementioned major design features, the *humeral component* also includes the following features to accept the *axle-pin* and *hearings A and B*: v-grooves at the base of the keyways to accept the *axle-pin*, a hole to accept *humeral bearing-A* and two keyways to accept *ulnar bearings-B*.




Image 11.5: Humeral Component, design features for mating to bearings and axle-pin

Ulnar Component

The *ulnar component* is made from wrought *Tivanium* alloy (Ti-6Al-4V) and has the following major design features: a *stem*, an *eye* and a *plasma spray region*.




Image 11.6: Ulnar Component, major design features



The *ulnar stem* is intended to be implanted with bone cement into the patient's ulnar medullary canal. It has a quadrangular cross-section, which is designed to prevent rotation of the *ulnar component* within the patient's medullary canal; and, it is curved along its proximal-distal axis to conform to the curved shaft of the proximal ulnar medullary canal.

The *ulnar eye* is both highly polished and nitrogen-enriched (as recommended in Section J, Titanium Articular Surface Coatings, of *Class II Special Controls Guidance Document: Knee Joint Patellofemorotibial and Femorotibial Metal/Polymer Porous-Coated Uncemented; Guidance for Industry and FDA) to limit wear of the apposing <i>Bearings A and B* upon which it articulates.

The proximal portion of the *ulnar stem* is covered with plasma spray, which is contained within the Coatings Characterization of Orthopedic Titanium 6Al4V Coatings Material Master File (FDA Ref: MAF-652, Letter of Authorization can be found in Section 22). The *plasma spray region* is designed to increase the coefficient of friction between the implant and the bone cement, thereby improving anti-rotational stability and aiding in resistance to movement of the *ulnar stem* with respect to the ulna.

The *ulnar component* is side-specific, meaning that there is a left *ulnar component* and a right *ulnar component* available for every size offering.

Humeral Bearing-A

The Humeral Bearing-A is made from Vivacit- E^{TM} (Vitamin E [α -tocopherol] stabilized, highly crosslinked ultra-high molecular weight polyethylene), which is contained within the Vitamin E Stabilized HXPE Material Master File authored by Zimmer (FDA Ref: MAF-1868). Humeral bearing-A is designed to articulate against the ulnar component throughout most of the humeroulnar joint range of motion, and has the following major design features: concave articulation surface, anterior and posterior ridges, a peg, four peripheral bumps and medial and lateral rails.




Image 11.7: Humeral Bearing-A, major design features



The humeral bearing-A concave articulation surface is intended to (along with the articulation surface of the ulnar bearings-B) allow for load sharing from the ulnar component onto humeral bearing-A throughout most of the range of motion of the humeroulnar joint. The posterior ridge specifically supports high-flexion motions (see Image 11.8, below). The humeral bearing-A anterior ridge is symmetrical to the posterior ridge in order to remove the human factors issue of having to discern anterior from posterior when assembling humeral bearing-A to the humeral component: humeral bearing-A cannot be put in backwards, because it is symmetrical.




Image 11.8: Posterior Ridge of Humeral Bearing-A supporting high flexion of the ulnar component with respect to the humeral component

The humeral bearing-A peg is designed to mate to the hole at the proximal base of the humeral yoke with a press (not a snap) fit, which is intended to limit motion between humeral bearing-A and the humeral yoke.

The *humeral bearing-A* (four) *peripheral bumps* are designed to mate against the interior walls of the *humeral yoke* with a press-fit, which is intended to limit motion (and therefore backside wear) of the *humeral bearing-A* relative to the *humeral yoke*.

The humeral bearing-A *medial and lateral rails* are designed to allow clearance for (as opposed to a direct mate with) *ulnar bearings-B*, which directly mate to the *humeral yoke* and to the *axle pin*.

Note: The *humeral bearing-A* may also be referred to as "humeral bearing" and "bearing A" throughout this submission.

Ulnar Bearings-B

The *ulnar bearings-B* (two per implant assembly) are made from Vivacit- E^{TM} (Vitamin E [a-tocopherol] stabilized, highly crosslinked ultra-high molecular weight polyethylene), which is contained within the Vitamin E Stabilized HXPE Material



Master File authored by Zimmer (FDA Ref: MAF-1868). Ulnar bearings-B are designed to articulate against the ulnar component throughout the entire humeroulnar joint range of motion, and have the following major design features: articulation surface, anti-rotation key, through-hole and compressible semi-circle.




Image 11.9: Ulnar Bearings-B, major design features

The ulnar bearings-B articulation surface is intended to (along with the articulation surface of the ulnar bearing-A) allow for load sharing from the ulnar component onto humeral bearing-B throughout the range of motion of the humeroulnar joint.

The ulnar bearings-B anti-rotation keys mate into the humeral yoke keyways with a press-fit to limit translational and rotational motion of ulnar bearings-B with respect to the humeral component.

The compressible semi-circle of one ulnar bearing-B is designed to align with the compressible semi-circle of the second ulnar bearing-B such that a nearly full circle is created. This design feature is intended to ensure proper alignment of ulnar bearings-B with respect to each other (during assembly of ulnar bearings-B to the axle pin using the ulnar bearing assembly tool) and, after proper assembly to the axle pin, to the humeral yoke keyways. Due to the fact that the compressible semi-circle is more easily deformed than the broad, apposing faces of ulnar bearings-B, it creates a press-fit between the ulnar bearings-B and between each ulnar bearing-B and the apposing interior walls of the humeral yoke. This press-fit is intended to limit mediolateral motion of the ulnar bearings-B along the axle-pin, thereby limiting backside wear of the ulnar bearings-B.

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017



Traditional 510(k) Premarket Notification

The *ulnar bearings-B through-hole* mates around the *axle pin* with a press-fit intended to limit motion of *ulnar bearings-B* with respect to the *axle pin*, and thus limit backside wear of *ulnar bearings-B*. The position of the *through-hole* is eccentric (with respect to the outer diameter of *ulnar bearings-B*) in order to place the largest cross-section of Vivacit-ETM material in opposition to the highest loads (Joint Reaction Forces, JRFs) induced throughout the flexion-extension cycle of the humeroulnar joint (see image 11.10, below).

Note: The *ulnar bearings-B* do NOT mate to or in any way interface with the *humeral screws*.

Note: The *ulnar bearings-B* may also be referred to as "ulnar bearing(s)" and "bearing(s) B" throughout this submission.




Image 11.10: Thickest Vitamin E HXPE placement to oppose highest JRFs

Axle Pin

The *axle pin* is made from *Zimaloy* (Cobalt-Chromium-Molybdenum Alloy), is designed to remain in a fixed position (not translate nor rotate) throughout the entire humeroulnar joint range of motion, and has the following major design feature: *v-groove*.




Image 11.11: Axle Pin, major design feature



The *v*-groove is designed to mate with the taper of the *humeral screws* that, upon screw thread engagement in the *humeral yoke*, rigidly compresses the *axle pin* against the *humeral yoke*.

Additionally, the *axle pin* is designed to rigidly mate with *ulnar bearings-B*, the *humeral screws* and the *humeral yoke*, creating a rigid assembly along which joint reaction forces are intended to be shared.

Note: The axle-pin may also be referred to as "axle" throughout this submission.

Humeral Screws

The *humeral screws* (two per implant assembly) are made from *Zimaloy* (Cobalt-Chromium-Molybdenum Alloy), are designed to remain in a fixed position throughout the entire humeroulnar joint range of motion, and have the following major design features: *headless*, *threads*, *pilot* and *taper*.




Image 11.12: Humeral Screws, major design feature

The humeral screw threads engage with threads of the humeral yoke while the humeral screw pilot engages a hole in the wall of the humeral yoke and the humeral screw taper engages with the axle pin v-groove, compressing the ends of the axle pin against the v-groove at the base of the humeral yoke keyway.

The *humeral screws* are *headless*, a feature that allows the *humeral screws* to countersink into the *humeral yoke*, which is intended to ensure the constant application of a clamp-load to the *axle-pin* and to avoid soft tissue impingement or irritation associated with a protruding screw head.

Note: While the *humeral screws* are designed to apply a consistent clamp-load to the *axle pin*, they are *not* directly loaded by any compressive joint reaction force (JRF) during the flexion-extension cycle (this design feature is intended to avoid screw back-out by not allowing any complete transverse load-reversal across the *humeral screw threads*). Instead, compressive JRFs are transmitted from the *ulnar implant*



component to the ulnar bearings-B to the axle pin to the v-grooves at the base of the humeral yoke keyways and finally throughout the humeral component. See Zimmer Research Report ZRR_WA 2598 12, New Zimmer Total Elbow Prosthesis Modular Connection Fatigue Testing for a more detailed explanation.

System Compatibility

No Zimmer Nexel Total Elbow component is compatible with any device outside the Zimmer Nexel Total Elbow system.

Size Interchangeability

All Zimmer Nexel Total Elbow humeral component sizes are designed to mate with all Zimmer Nexel Total Elbow ulnar component sizes. The axle pin size 4, humeral bearing-A size 4, and ulnar bearings-B size 4 are designed to mate with all size 4 humeral components. The axle pin size 5/6, humeral bearing-A size 5/6, and ulnar bearings-B size 5/6 are designed to mate with all size 5 and size 6 humeral components. The humeral screws are designed to mate with all Zimmer Nexel Total Elbow humeral components and axle pin sizes.

Engineering Drawings

Representative <u>engineering drawings</u> for the proposed device are included at the end of this section.



. s.

Catalog Numbers

All catalog numbers for the Zimmer Nexel Total Elbow are listed below.

Catalog Number	Description
00-8400-014-07	Zimmer Nexel Total Elbow - Ulnar Stern, Size 4, 75 mm, Left
00-8400-014-11	Zimmer Nexel Total Elbow - Ulnar Stem, Size 4, 115 mm, Left
00-8400-024-07	Zimmer Nexel Total Elbow - Ulnar Stem, Size 4, 75 mm, Right
00-8400-024-11	Zimmer Nexel Total Elbow - Ulnar Stem, Size 4, 115 mm, Right
00-8400-015-07	Zimmer Nexel Total Elbow - Ulnar Stem, Size 5, 75 mm, Left
00-8400-015-11	Zimmer Nexel Total Elbow - Ulnar Stem, Size 5, 115 mm, Left
00-8400-025-07	Zimmer Nexel Total Elbow - Ulnar Stem, Size 5, 75 mm, Right
00-8400-025-11	Zimmer Nexel Total Elbow - Ulnar Stem, Size 5, 115 mm, Right
00-8400-016-09	Zimmer Nexel Total Elbow - Ulnar Stem, Size 6, 90 mm, Left
00-8400-016-11	Zimmer Nexel Total Elbow - Ulnar Stem, Size 6, 115 mm, Left
00-8400-026-09	Zimmer Nexel Total Elbow - Ulnar Stem, Size 6, 90 mm, Right
00-8400-026-11	Zimmer Nexel Total Elbow - Ulnar Stem, Size 6, 115 mm, Right
00-8400-044-10	Zimmer Nexel Total Elbow - Humeral Stern, Size 4, 100 mm
00-8400-044-15	Zimmer Nexel Total Elbow - Humeral Stem, Size 4, 150 mm
00-8400-045-10	Zimmer Nexel Total Elbow - Humeral Stem, Size 5, 100 mm
00-8400-045-15	Zimmer Nexel Total Elbow - Humeral Stem, Size 5, 150 mm
00-8400-055-15	Zimmer Nexel Total Elbow - Humeral Stem, Long Flange, Size 5, 150 mm
00-8400-055-20	Zimmer Nexel Total Elbow - Humeral Stem, Long Flange, Size 5, 200 mm
00-8400-046-10	Zimmer Nexel Total Elbow - Humeral Stem, Size 6, 100 mm
00-8400-046-15	Zimmer Nexel Total Elbow - Humeral Stern, Size 6, 150 mm
00-8400-056-15	Zimmer Nexel Total Elbow - Humeral Stem, Long Flange, Size 6, 150 mm
00-8400-056-20	Zimmer Nexel Total Elbow - Humeral Stem, Long Flange, Size 6, 200 mm
Kit # 00-8400-094-00 00-8400-194-01 00-8400-194-02	Zimmer Nexel Total Elbow - Articulation Kit, Size 4 Kit Contents: 1 Humeral Bearing-A Size 4 2 Ulnar Bearings- B Size 4
00-8400-194-04	1 Axle Pin Size 4
Kit # 00-8400-095-00 00-8400-195-01 00-8400-195-02 00-8400-195-04	Zimmer Nexel Total Elbow - Articulation Kit, Size 5/6 Kit Contents: 1 Humeral Bearing-A Size 5/6 2 Ulnar Bearings- B Size 5/6 1 Axle Pin Size 5/6
Kit # 00-8400-090-00	Zimmer Nexel Total Elbow - Humeral Screw Kit Kit Contents:
00-8400-190-00	2 Humeral Screws

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

2

.



Materials

The humeral and ulnar components are made of wrought Tivanium alloy (Ti-6AI-4V) and are titanium plasma sprayed (on regions specified in the engineering drawings at the end of this section). The bearings A and B are made of Vitamin E Highly Crosslinked Polyethylene). The axle pin and humeral screws are made of Zimaloy (Cobalt-Chromium-Molybdenum Alloy).

The following tests contained within the Material Master File for Vitamin E Stabilized Highly Crosslinked Polyethylene (Master File MAF-1868) are applicable to usage in elbow devices and support substantial equivalence of the Vitamin E Stabilized HXPE material when used in elbow devices:

Report Number	Report Name	Master File Location
ZRR_WA_2401_11 Rev. 2	TESTING OF VITAMIN E HIGHLY CROSSLINKED POLYETHYLENE PROPERTIES PER ASTM F2759-09	Amendment 02, Exhibit B
N102082	NAMSA - TOXICOLOGICAL RISK ASSESSMENT	Volume 1, Page 54
T0118_913/S	13 WEEK SYSTEMIC TOXICITY STUDY IN RATS FOLLOWING SUBCUTANEOUS IMPLANTATION	Volume 1, Page 86
V0014_130	CYTOTOXICITY STUDY USING THE ISO ELUTION METHOD	Volume I, Page 76
V0023_211	BACTERIAL REVERSE MUTATION STUDY - EXTRACTS	Volume 1, Page 108
V0573_000/S	GENOTOXICITY: MOUSE LYMPHOMA ASSAY	Volume 1, Page 125
T0566_500, T0566_501	MOUSE PERIPHERIAL BLOOD MICRONUCLEUS STUDY	Volume 1, Pae 147
T1251_800	ISO INTRACUTANEOUS STUDY IN RABBITS	Volume 1, Page 179
T1261_300	ISO GUINEA PIG MAXIMIZATION SENSITIZATION TEST	Volume I, Page 189
T1250_812	ISO MUSCLE IMPLANTATION STUDY IN RABBITS - 12 WEEKS	Volume 1, Page 205
T1250_802	ISO MUSCLE IMPLANTATION STUDY IN RABBITS - 2 WEEKS	Volume 1, Page 217
ZRR_WA_2338_11	VITAMIN E HIGHLY CROSSLINKED POLYETHYLENE SHELF LIFE STUDY	Volume I, Page 230
T02625-500	ISO SYSTEMIC TOXICITY STUDY IN MICE	Volume 1, Page 245
T0118_926	26 WEEK SYSTEMIC TOXICITY STUDY IN RATS FOLLOWING SUBCUTANEOUS IMPLANTATION	Volume I, Page 256
ZRR_WA_2424_11	EVALUATION OF VITAMIN E DOPED AND CROSSLINKED UHMWPE WEAR DEBRIS IN THE RABBIT KNEE FOLLOWING PERCUTANEOUS INJECTION	Amendment 01, Page 9



Report Number	Report Name	Master File Location	
10T_60283_02	BIOLOGICAL RESPONSE TO WEAR DEBRIS Note: The particulate study was done to detect local and systemic response to vitamin E grafted UHMWPE. A rabbit model which is widely accepted was used. The particulates were produced from production materials and have size distribution that is representative of those that have been isolated from tissues harvested during hip or knee revision surgeries. The particles were injected directly into the knee capsule of the animals. This study was not focused on any particular implant or device. The emphasis was on inflammatory response to the particulates. The study covered all devices that generate particulates.	Volume 2, Page 243	
ZRR_WA_2409_11	VIATMIN E STABILIZED HIGHLY CROSSLINKED POLYETHYLENE EXTENDED AGING STUDY	Volume 2, Page 394	
ZRR_WA_2373_11	R_WA_2373_11 ENVIRONMENTAL STRESS CRACKING OF VITAMIN E, LONGEVITY, CONVENTIONAL UHMWPE		
ZRR_WA_2403_11	URR_WA_2403_11 VITAMIN E STABILIZED HIGHLY CROSSLINKED POLYETHYLENE- EXTRACTABILITY STUDY		
ZRR_WA_2412_11	VITAMIN E STABILIZED HIGHLY CROSSLINKED POLYETHYLENE- PRESSURE INDUCED LEACHING	Volume 2, Page 492	
ZRR_WA_2382_11	VITAMIN E STABILIZED HIGHLY CROSSLINKED POLYETHYLENE- TRANSFORMATION PRODUCTS EVALUATION BY HPLC AND GC	Volume 2, Page 507	
ZRM_WA_0293_J1	LAMELLAE THICKNESS	Volume 2, Page 566	
ZRM_WA_0277_11	MOLECULAR WEIGHT	Volume 2, Page 578	
CPG REPORT_ 11622_1_	GC/MS ANALYSIS OF UHMWPE	Amendment 02, Exhibit A	

Surface Characteristics

The eye of the ulnar component is treated with a nitrogen ion implantation (also known as nitriding) process (to increase *Tivanium* alloy surface hardness) on the region specified in the <u>engineering drawings</u> at the end of this section.

The articulating surfaces of the *ulnar component* are polished to an R_a value of ((b)(4)) or less and the articulating surfaces of the *bearings A and B* are machined to an R_a value of (b)(4) or less to provide for smooth articulating surfaces.



Mechanism of Action

The Zimmer Nexel Total Elbow implant assembly is a constrained prosthesis designed to allow anatomic reconstruction of the humeroulnar joint with the natural change in humeroulnar joint center of rotation throughout the flexion/extension arc based on research presented in the Zimmer Research Report ZRR_WA_2476_11, *Biomechanics Rationale for the Loading of Total Elbow Prosthesis During In Vitro Simulations*. See table below for range of motion details.

		Feature	Proposed Device: The Zimmer Nexel Total Elbow
ent)		Varus/Valgus Angular Laxity	Neutral to Maximum Varus: Neutral to Maximum Valgus: Total Included Angle:
on is of ovement Component)		Internal/External Rotational Laxity	Neutral to Maximum Internal Rotation: (b) Neutral to Maximum External Rotation: (4) Total Included Angle: (b)
of Motion in terms onent mov umeral C	Proximal/Distal Translational Laxity Proximal/Distal Translational Laxity Anterior/Posterior Translation Laxit		Neutral to Maximum Proximal: (b)(4) Neutral to Maximum Distal: Total Proximal/Distal Laxity:
			Neutral to Maximum Anterior: Neutral to Maximum Posterior: Total Anterior/Posterior Laxity:
Ra (del Ulnar C with respect		Medial/Lateral Translation Laxity	Neutral to Maximum Medial: Neutral to Maximum Lateral: Total Medial/Lateral Laxity:
×		Flexion/Extension Range (defined in terms of degrees flexion)	All Humeral and Ulnar size combinations: (b)(4)

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118



Surgical Instrumentation Unique to the Device

Following are the instruments unique to the Zimmer Nexel Total Elbow. These are Class II accessories to the device. Applicable <u>material standards</u> are listed in a table at the end of this section.

	The Zimmer Nexel Total Elbow – Class II Accessories to the Implant (fable Continued.on Next Page)								
Cataloğ Number	Description	Material	Tissue Contact?	Imáge (nót (ö]scale)	Întended [®] Use				
00-8401- 002-00	T-Handle	17-4 SS 303 SS 420 SS 316 SS	yes	Comparison of the second se	Quick-connect handle for use with humeral and ulnar rasps.				
00-8401- 008-00	Reciprocating Saw Adapter	303 SS 420 SS 316 SS	yes	Contraction of the second s	Adapter to connect humeral and ulnar rasps to a compatible large bone reciprocating saw via the single use connector.				
00-8401- 008-01	Single Use Connector	455 SS	yes, skin only		Connects the reciprocating saw adaptor to a compatible large-bone reciprocating saw.				
00-8401- 009-00	Slide Hammer	17-4 SS 304 SS 17-4 SS	yes		Mates with humeral and ulnar stem extractors to deliver impact blows when extractor component is attached to the implant.				
00-8401- 018-00	Ulnar Bearing Tamp	17-4 SS	yes		To compress the articulating assembly (ulnar bearings and axle pin) into the humeral implant. (An alternate instrument to the articulation inserter.)				
00-8401- 012-04	Trephine Stabilizer, Size 4	17-4 SS	yes	P	Guides the finishing cut of the trephine through the				
00-8401- 012-05	Trephine Stabilizer, Size 5	17-4 SS	yes	anterior aspec	anterior aspect of the humerus.				

.



	The Zimmer Nexel Total Elbow – Class II Accessories to the Implant (Table Continued from Previous Page, Continued on Next Page)							
Catalog Number	Description	Material	Tissue Contact?	Image (not to scale)	Intended Use			
00-8401- 014-07	Ulnar Provisional, Size 4, 75 mm, Left	17-4 SS	yes					
00-8401- 014-11	Ulnar Provisional, Size 4, 115 mm, Left	17-4 SS	yes					
00-8401- 015-07	Ulnar Provisional, Size 5, 75 mm, Left	17-4 SS	yes					
00-8401- 015-11	Ulnar Provisional, Size 5, 115 mm, Left	17-4 SS	yes					
00-8401- 016-09	Ulnar Provisional, Size 6, 90 mm, Left	17-4 SS	yes					
00-8401- 016-11	Ulnar Provisional, Size 6, 115 mm, Left	17-4 SS	yes		preparation of the	To assess the bone preparation of the ulnar		
00-8401- 024-07	Ulnar Provisional, Size 4, 75 mm, Right	17-4 SS	yes		medullary canal, and to allow an intraoperative range of motion.			
00-8401- 024-11	Ulnar Provisional, Size 4, 115 mm, Right	17-4 SS	yes					
00-8401- 025-07	Ulnar Provisional, Size 5, 75 mm, Right	17-4 SS	yes					
00-8401- 025-11	Ulnar Provisional, Size 5, 115 mm, Right	17-4 SS	yes					
00-8401- 026-09	Ulnar Provisional, Size 6, 90 mm, Right	17-4 SS	yes					
00-8401- 026-11	Ulnar Provisional, Size 6, 115 mm, Right	17-4 SS	yes					
00-8401- 019-00	Articulation Inserter	17-4 SS	yes		To compress the articulating assembly (ulnar bearings and axle pin) and ulnar implant into the humeral implant. Also, to compress the humeral bearing into the humeral implant. (This is an alternate instrument to the ulnar bearing tamp and the humeral bearing driver.)			
00-8401- 028-00	Ulnar Stem Inserter	Radel	yes		To transmit force to insert the ulnar implant without damaging the articular surface.			

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118



.

Т	The Zimmer Nexel Total Elbow – Class II Accessories to the Implant (Table Continued from Previous Page, Continued on Next Page)						
Catalog Number	Description	Material	Lissue Contact?	Intage (not to scale)	Intended Use		
00-8401- 029-00	Implant Extractor Hook	17-4 SS	yes		To transmit force from the slide hammer to an implanted ulnar implant or humeral implant to facilitate their removal. The extractor hook is attached to the slide hammer by a threaded connection.		
00-8401- 033-01	Ulnar Rasp, Pilot, Left	17-4 SS	· yes		Removes bone from		
00-8401- 033-02	Ulnar Rasp, Pilot, Right	17-4 SS	yes		the ulnar medullary canal and olecranon and shapes the medullary canal to		
00-8401- 034-01	Ulnar Rasp, Size 4/5, Left	17-4 SS	yes	4	create a precise geometry for a size- matched ulnar implant.		
00-8401- 034-02	Ulnar Rasp, Size 4/5, Right	17-4 SS	yes		Can be mated to the reciprocating saw adaptor or to the T-		
00-8401- 036-01	Ulnar Rasp, Size 6, Left	17-4 SS	yes		handle.		
00-8401- 036-02	Ulnar Rasp, Siz e 6, Right	17-4 SS	yes				
00-8401- 039-00	Ulnar Clearance Template	17-4 SS	ycs		To check for adequate bone removal on the medial and lateral sides of the ulna for clearance for the articulation insertion instruments.		



:	The Zinimer Nexel Total Elbow – Class II Accessories to the Implant (Table Continued from Previous Page, Continued on Next Page)							
Catalog Number	Description	Material	Tissue Contact?	Image (not to scale)	Intended Use			
00-8401- 044-10	Humeral Provisional, Size 4, 100 mm	17-4 SS	yes					
00-8401- 044-15	Humeral Provisional, Size 4, 150 mm	17-4 SS	yes					
00-8401- 045-10	Humeral Provisional, Size 5, 100 mm	17-4 SS	yes					
00-8401- 045-15	Humeral Provisional, Size 5, 150 mm	17-4 SS	yes		Dimensionally			
00-8401- 046-10	Humeral Provisional, Size 6, 100 mm	17-4 SS	yes	BE	represents the humeral implant to help assess the preparation of the			
00-8401- 046-15	Humeral Provisional, Size 6, 150 mm	17-4 SS	yes		humeral medullary canal, and to intraoperatively check range of motion.			
00-8401- 055-15	Humeral Provisional, Long Flange, Size 5, 150mm	17-4 SS	yes					
00-8401- 055-20	Humeral Provisional, Long Flange, Size 5, 200mm	17-4 SS	yes					
00-8401- 056-15	Humeral Provisional, Long Flange, Size 6, 150mm	17-4 SS	yes					
00-8401- 056-20	Humeral Provisional, Long Flange, Size 6, 200mm	17-4 SS	yes					
00-8401- 058-04	Humeral Stem Inserter, Size 4	17-4 SS	yes	A	To insert the humeral implant into the humeral canal without			
00-8401- 058-05	Humeral Stem Inserter, Size 5/6	17-4 SS	yes	PO	humeral canal without damaging the humeral bearing.			
00-8401- 059-04	Humeral Extractor Plate, Size 4	17-4 SS	yes		To transmit force from slide hammer to humeral implant in order to facilitate			
00-8401- 059-05	Humeral Extractor Plate, Size 5/6	17-4 SS	yes		removal. Extractor plate is attached to humeral implant by the humeral extractor screws.			



	The Zimmer Nexel Total Elbow – Class II Accessories to the Implant (Table Continued from Previous Page, Continued on Next Page)							
Catalog Number	Description	Material	Tissue Contact?	Image (not to scale).	Intended Use			
00-8401- 060-00	Humeral Awl Reamer	17-4 SS	yes	×°	To open the humeral medullary canal. Also, provides a means of ensuring the trochlear resection is adequate to accept subsequent humeral rasps.			
00-8401- 061-00	I/E Alignment Rod	17-4 SS	yes		To assess and adjust, if necessary, the internal/external rotation of the humeral and ulnar rasps.			
00-8401- 064-00	Humeral Rasp, Pilot	17-4 SS	yes		Removes bone from th c			
00-8401- 064-10	Humeral Rasp, Size 4, 100 mm	17-4 SS	yes		humeral medullary canal to create a precise geometry			
00-8401- 064-15	Humeral Rasp, Size 4, 150 mm	17-4 SS	yes		for a size-matched humeral implant. Humeral rasps also used to orient and securely			
00-8401- 065-10	Humeral Rasp, Size 5, 100 mm	17-4 SS	yes	P	support the humeral cut guide, and to guide the			
00-8401- 065-15	Humeral Rasp, Size 5, 150 mm	17-4 SS	yes		trephine cut.			
00-8401- 066-10	Humeral Rasp, Size 6, 100 mm	17-4 SS	yes					
00-8401- 066-15	Humeral Rasp, Size 6, 150 mm	17-4 SS	yes		· · ·			



	The Zimmer Nexel Total Elbow – Class H Accessories to the Implant (Table Continued from Previous Page, Continued on Next Page)							
Catalog Number	Description	Material	Tissue Contact?	Image (not to scale)	Intended Use			
00-8401- 070-45	Flexible Solid Reamer, 4.5 mm	1 3-8 SS	yes	1	To remove bone from the ulnar medullary canal. Connects to a			
00-8401- 070-05	Flexible Solid Reamer, 5 mm	13-8 SS	yes		surgical drill or reamer driver.			
00- 840 1- 072-01	Flexible Cannulated Reamer, 5.5 mm	13-8 SS	yes					
00-8401- 072-02	Flexible Cannulated Reamer, 6 mm	13-8 SS	yes					
00-8401- 072-03	Flexible Cannulated Reamer, 6.5 mm	1 3-8 S S	yes	A A				
00-8401- 072-04	Flexible Cannulated Reamer, 7 mm	13-8 SS	yes	to remove bone from t and humeral medullar Connects to a surgical	Mates over the guide wire. Used			
00-8401- 072-05	Flexible Cannulated Reamer, 8 mm	13-8 SS	yes		to remove bone from the ulnar and humeral medullary canals. Connects to a surgical drill or			
00-8401- 072-06	Flexible Cannulated Reamer, 9 mm	13-8 SS	yes			reamer driver.		
00-8401- 072-07	Flexible Cannulated Reamer, 10 mm	13-8 SS	yes					
00-8401- 072-08	Flexible Cannulated Reamer, 11 mm	13-8 SS	yes					
00-8401- 072-09	Flexible Cannulated Reamer, 12 mm	1 3-8 SS	yes					
00-8401- 074-00	Trephine, Size 4	17-4 SS 13-8 SS	yes		Removes bone from the distal humerus to accept a given size-			
00-8401- 075-00	Trephine, Size 5/6	17-4 SS 13-8 SS	yes		matched humeral implant. Connects to a surgical drill.			
00-8401- 078-00	Humeral Bearing Driver	17-4 SS	yes		To deliver a compressive force sufficient to fully seat the humeral bearing into the humeral implant when used with the humeral bearing driver pin. Also used to install a humeral bearing component into the humeral implant during a revision surgery. (This is an alternate instrument to the articulation inserter.)			



The Zimmer Nexel Total Elbow – Class II Accessories to the Implant (Table Continued from Previous Page, Continued on Next Page)							
Catalog Number	Description	Material	Tissue Contact?	Image (not to scale)	Intended Use		
00- 84 01- 079-00	Humeral Bearing Driver Pin	440A SS	yes		Subcomponent of the humeral bearing driver. Can also be used to aid in internal/external alignment of the ulnar provisionals.		
00-8401- 081-00	Ulnar Bearing Assembly Tool	17-4 SS 455 SS 316 SS	yes	R	To assemble the axle pin to ulnar bearings-B; and, to assemble the ulnar bearings-B onto the axle pin through the proximal head of the ulnar Implant. Performed by compressing the handle of the tool manually. To be used in situ.		
00-8401- 082-00	Humeral Bearing Placement Tool	420 SS 304 SS	yes	Å	To place humeral bearing-A into the base of the humeral implant yoke.		
00-8401- 084-00	Humeral Screw Holder, Size 4	17-4 SS	yes		To securely hold the humeral screws and align the screws to a size-matched		
00-8401- 085-00	Humeral Screw Holder, Size 5/6	17-4 SS	yes		humeral implant.		
00-8401- 088-00	Humeral Screw Installation Driver Bit	440C SS LDPE	yes		Connects to the humeral screw installation driver to install the humeral screws into the humeral implant, and prevents the user from over-torquing the humeral screws. An alternative to the torque-limiting humeral screw driver, below. NOTE: Cannot mate with the humeral screw removal driver.		
00-8401- 089-00	Humeral Screw Removal Driver	455 SS 303 SS Silicone	yes		To remove the humeral screws from the humeral implant only when mated with the humeral screw removal driver bit. NOTE: Will <i>not</i> allow installation torque to be applied.		
00-8401- 091-00	Humeral Bearing Extractor	17-4 SS	yes	Q	To remove the humeral bearing from the humeral implant.		
00-8401- 092-00	Articulation Extractor	17-4 SS	yes		To remove the articulation assembly (ulnar bearings and axle pin) and ulnar implant from the humeral implant.		

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118



The Zimmer Nexel Total Elbow – Class II Accessories to the Implant (Table Continued from Previous Page, End of Table)										
Catalog Number	Description	Material	Tissue Contact?	lmage (not to scale)	Intended lise					
00-8401- 093-00	Humeral Extractor Screw	455 SS	yes	Contract of the second	Attaches the humeral extractor plate to the humeral implant to facilitate humeral component extraction.					
00-8401- 097-00	Humeral Screw Removal Driver Bit	440C SS 316 SS	yes		Used with the humeral screw removal driver. NOTE: Cannot mate with the humeral screw installation driver.					
00-8401- 004-00	Humeral Cut Guide, Size 4	17-4 SS	yes		Connects to humeral rasps. Used to facilitate the creation of parallel cuts equidistant from the longitudinal axis of the rasp when used with a surgical blade.					
00-8401- 005-00	Humeral Cut Guide, Size 5/6	17-4 SS	yes	U.						
00-8401- 080-00	Torque-Limiting Humeral Screw Driver	440C SST Ultem	yes		To install and remove the humeral screws into and from the humeral implant. Prevents the user from over- torquing the humeral screws upon insertion. This instrument is intended to be used for a single surgery only.					

;



Surgical Instrumentation Not Unique to the Device

Following are the instruments not unique to the Zimmer Nexel Total Elbow. These are Class I instruments currently able to be used with other surgical devices cleared for market in the United States.

The Zimmer Nexel Total Elbow – Class I Instruments				
Catalog Number	Description			
00-4812-035-00	Small Hex Screwdriver			
47-2255-008-00	Sterile Ball Tip Guide Wire 2.4 x 70			
31-8106-168-00	Ulnar Awl Reamer			
00-4811-035-01	Humeral Screw Installation Driver			
00-5049-053-00	Quick Use Curette			



•

Traditional 510(k) Premarket Notification

Instrument Material Standards

The Zimme	r Nexel Total Elbow: Instrument Material Standards				
Material	Applicable Standards per Zimmer Engineering Specificati				
6061-T6 Aluminum	ASTM B209, ASTM B211, ASTM B221				
13-8 SS	ASTM A484, ASTM A555, ASTM A564				
17-4 SS	ASTM A484, ASTM F899				
18-8 SS	ASTM A484, ASTM A555, ASTM A580				
303 SS	ASTM F899				
304 SS	ASTM A276, ASTM A484				
316 SS	ASTM F899				
316L SS	ASTM F138, ASTM F139, ASTM 1350, ISO 5832-1				
420 SS	ASTM F899				
440A SS	ASTM F899, ASTM A555				
440C SS	ASTM F899				
455 SS	ASTM F899, ISO 7153-1				
Phenolic	No Standard Applicable				
L 605 (Co-Cr-W-Ni alloy)	ASTM F90				
LDPE	ASTM D1248				
Radel	No Standard Applicable				
Silicone	No Standard Applicable				
Ultem	ASTM D5205				




.

Traditional 510(k) Premarket Notification

Methods, Facilities and Controls

Method of Manufacturing





Packaging

Humeral Stem



The above listed materials are identical to those currently in use for Zimmer sterile implant packages with the only exception being that the configuration of the interior components correspond to the particular product being packaged.



The above listed materials are identical to those currently in use for Zimmer sterile implant packages with the only exception being that the configuration of the interior components correspond to the particular product being packaged.



Articulation Kit



The above listed materials are identical to those currently in use for Zimmer sterile implant packages with the only exception being that the configuration of the interior components correspond to the particular product being packaged.

Humeral Screw Kit



The above listed materials are identical to those currently in use for Zimmer sterile implant packages. Only the configuration of the interior components changes to correspond to the particular product being packaged.

Pyrogenicity

This device is not labeled as nonpyrogenic. Per USP XXIII (161), requirements for specified endotoxin levels do not apply to orthopaedic implants.

Latex

There is no natural latex rubber in this product or its packaging.

Records processed under FOIA request #20/16 4653 Released by CDRH on 07/16/2017 2/27/13 email



K123862 Zimmer Nexel Total Elbow Refuse to Accept – Zimmer Response

FDA Question 1

			Yes	N/A	No	
4.	Eit	bmission contains 510(k) Summary or 510(k) Statement her a) or b) must be answered "Yes" to be considered nplete. Identify any missing element(s) as Comments.	X			
	a.	Summary contains all elements per 21 CFR 807.92 See also 510(k) Summary Checklist			x	
	b.	Statement contains all elements per 21 CFR 807.93				
	Comments: The 510(k) summary is missing a complete Device Description, (lacking explanation of how the device functions, the scientific concepts that form the basis for the device, and the significant physical and performance characteristics of the device, such as device design, material used, and physical properties) and is missing a comparison to predicate devices.					

2

, i



Zimmer Response to FDA Question 1

Please see Attachment 1 for a replacement 510(k) summary.

Within the device description section of the replacement 510(k) summary, the following verbiage was added to that of the original submission:

How the Device Functions: The Zimmer Nexel Total Elbow is an implant designed to replace the articulating surfaces of and restore motion to the humeroulnar (elbow) joint. The implant is a constrained device assembly and consists of the following components: humeral component, ulnar component, humeral bearing-A, ulnar bearings-B, axle pin and humeral screws.

The Following Scientific Concepts, Design Features and Physical Properties form the Basis for the Zimmer Nexel Total Elbow: The humeral component has a humeral stem designed to be implanted with bone cement into the patient's humeral medullary canal, an anterior flange designed to accept a bone graft and limit torsional and posterior migration, a humeral yoke with rounded corners to avoid the creation of stress risers within the medial and lateral humeral supracondylar columns, and plasma spray region to enhance fixation to bone cement within the medullary canal, and to improve fatigue strength. The ulnar component has an ulnar stem designed to be implanted with bone cement into the patient's ulnar medullary canal, an ulnar eye that is both highly polished and nitrogen-enriched to limit wear of the apposing polymer bearings, and plasma spray region to enhance fixation to bone cement within the medullary canal. Bearings A and B are designed to broadly distribute joint reaction forces.

Materials Used: The humeral and ulnar components are made of a titanium alloy, the bearings A and B are made of Vitamin E highly cross-linked polyethylene, and the humeral screws and axle pin are made of a cobalt-chromium-molybdenum alloy.

Comparison to the Predicate: The proposed device (Zimmer Nexel Total Elbow) and the predicate device (Coonrad/Morrey Total Elbow) have the same intended use and similar indications for use. The proposed device humeral and ulnar components are very similar in terms of materials used and design/dimensions to the predicate device humeral and ulnar components. The primary differences between the proposed and predicate devices are that the proposed device bearing components are made of Vitamin E highly cross-linked polyethylene (the predicate device bearing components are made of ultrahigh molecular-weight polyethylene), and the proposed device has bearings that articulate on both the outer and inner diameters of the ulnar eye (the predicate device has bearings that articulate on the inner diameter of the ulnar eye). Additionally, the proposed device uses screws and a 1-part axle pin to fix the bearings in place, while the predicate device uses a 2-part (snap-fit) axle pin.



FDA Question 2

		Yes	N/A	No
17.	Submission includes an analysis of why any differences between the subject device and predicate(s) do not render the device NSE (e.g., do not constitute a new intended use, affect safety or effectiveness, or raise different questions of safety and effectiveness) (see section 513(i)(1)(A) of the FD&C Act)			X
	If there is no difference between the subject and predicate(s) with respect to indications for use or technology, this should be explicitly stated, in which case "N/A" should be selected. Select "No" only if the submission does not include an analysis of differences as described above or a statement that there are no differences. Note that due to potential differences in manufacturing that may not be known to the submitter, no identified differences does not necessarily mean that no performance testing is needed.			
	Comments: The sponsor provides a tabulated comparison of and predicate devices. From this tabulated list, it is clear the differences between the subject and predicate devices with indications for use and design features. The sponsor does no analysis/discussion of why these differences do not render the	ere ar espec of prov	e t to vide aı	1



Zimmer Response to FDA Question 2

Please see Attachment 2 for a replacement Substantial Equivalence Discussion.

Within the replacement Substantial Equivalence Discussion, the following verbiage was added to address FDA Question 2:

Why Any Differences Do Not Render the Device NSE

The major technological differences between the proposed device (Zimmer Nexel Total Elbow) and the predicate device (the Zimmer Coonrad/Morrey Total Elbow) are listed below along with the performance tests conducted to demonstrate substantial equivalence (the listed performance tests and their conclusions are described in more detail in *Section 18a Performance Testing – Bench* of this premarket notification submission):

Device Component	Major Difference Between Proposed and Predicate Devices	Performance Test Conducted to Demonstrate Substantial Equivalence
Humeral Component	 Distal design features for implant assembly Plasma Spray 	 Durability Testing (ZRR_WA_2542_12) Modular Connection Fatigue Testing (ZRR_WA_2598_12) Fatigue Testing (ZRR_WA_2554_12 & ZRR_WA_2589_12)
Ulnar Component	- Nitrogen enrichment	 Wear Testing (ZRR_WA_2407_11 & ZRR_WA_2552_12) Fatigue Testing (ZRR_WA_2585_12, ZRR_WA_2615_12, ZRR_WA_2615_12, ZRR_WA_2554_12 & ZRR_WA_2626_12)
Polyethylene Components	 Polyethylene material Surfaces against which components are intended to articulate 	 Wear Testing (ZRR_WA_2407_11 & ZRR_WA_2552_12) Durability Testing (ZRR_WA_2542_12)
Axle Pin	- Single Solid vs. 2-Part Hollow construction	 Wear Testing (ZRR_WA_2407_11 & ZRR_WA_2552_12) Durability Testing (ZRR_WA_2542_12) Modular Connection Fatigue Testing (ZRR_WA_2598_12)
Humeral Screws	- Existence of screws	 Wear Testing (ZRR_WA_2407_11 & ZRR_WA_2552_12) Durability Testing (ZRR_WA_2542_12) Modular Connection Fatigue Testing (ZRR_WA_2598_12)



Current questions of safety and effectiveness (clinically observed failure modes) associated with elbow prostheses are humeral and ulnar component loosening, humeral and ulnar component fracture, polyethylene wear, system (implant assembly) durability and system (implant assembly) disassociation (also known as modular connection fatigue). To demonstrate that the proposed device is at least as safe and effective as the predicate device despite the compounding effects of these technological differences and previously outlined dimensional differences, fatigue (to assess resistance to fracture), wear, durability and modular connection fatigue performance tests were conducted (testing conclusions are described in more detail in *Section 18a Performance Testing – Bench* of this premarket notification submission). Stem loosening was not the focus of a performance test for the New Zimmer Total Elbow, because the variables (cementing technique, stem geometry, stem material and stem surface finish) associated with stem loosening are consistent from the predicate to the proposed device; and, therefore the proposed device (the New Zimmer Total Elbow) is expected to be as safe and effective as the predicate device (the Coonrad/Morrey Total Elbow) in terms of ability to resist stem loosening.

As reflected in the dFMEA documentation associated with the proposed device, the Zimmer Nexel Total Elbow does not raise any new/different types of questions of safety or effectiveness.

The proposed and predicate devices have the exact same intended use. The proposed and predicate devices share similar indications for use; where, the only changes in wording of the indications for use are that the proposed device has indications for use that are intended to add specificity to the indications for use of the proposed device. This added specificity does not change the intended meaning of the indications for use of the proposed device as compared to those of the predicate.

Because the aforementioned differences between the proposed and predicate devices do not constitute a new intended use, affect safety or effectiveness or raise different questions of safety and effectiveness, these differences *do not* render the proposed device NSE (not substantially equivalent).



FDA Question 3

			Yes	N/A	No
27.	Assessment of the need for sterilization information				
	a.	Identification of device, and/or accessories, and/or components that are provided sterile.	X		
	b.	Identification of device, and/or accessories, and/or components that are end user sterilized			Х
	ins pro re- coi	mments: The sponsor has identified a number of patient c truments; however, the sponsor has not discussed whethe ovided sterile or non-sterile and whether they are able to b sterilized. Therefore, please provide 1) a discussion as to nponents are end user sterilized and 2) instructions for cl rilization for these components.	r they : be clear which	are ned an	



Zimmer Response to FDA Question 3

Of the Class II Accessories to the Implant listed in Section 11a. (Device Description) of the premarket notification submission K123862, the following are provided non-sterile and are reusable and end-user sterilized (Table 1, below). Instructions for cleaning and re-sterilization of these instruments are provided in <u>Attachment 3</u>.

Catalog Number	Description	
00-8401-002-00	T-Handle	
00-8401-008-00	Reciprocating Saw Adapter	
00-8401-008-01	Single Use Connector	
00-8401-009-00	Slide Hammer	
00-8401-018-00	Ulnar Bearing Tamp	
00-8401-012-04	Trephine Stabilizer, Size 4	
00-8401-012-05	Trephine Stabilizer, Size 5	
00-8401-014-07	Ulnar Provisional, Size 4, 75 mm, Left	
00-8401-014-11	Ulnar Provisional, Size 4, 115 mm, Left	
00-8401-015-07	Ulnar Provisional, Size 5, 75 mm, Left	
00-8401-015-11	Ulnar Provisional, Size 5, 115 mm, Left	
00-8401-016-09	Ulnar Provisional, Size 6, 90 mm, Left	
00-8401-016-11	Ulnar Provisional, Size 6, 115 mm, Left	
00-8401-024-07	Ulnar Provisional, Size 4, 75 mm, Right	
00-8401-024-11	Ulnar Provisional, Size 4, 115 mm, Right	
00-8401-025-07	Ulnar Provisional, Size 5, 75 mm, Right	
00-8401-025-11	Ulnar Provisional, Size 5, 115 mm, Right	
00-8401-026-09	Ulnar Provisional, Size 6, 90 mm, Right	
00-8401-026-11	Ulnar Provisional, Size 6, 115 mm, Right	
00-8401-019-00	Articulation Inserter	
00-8401-028-00	Ulnar Stem Inserter	
00-8401-029-00	Implant Extractor Hook	
00-8401-033-01	Ulnar Rasp, Pilot, Left	
00-8401-033-02	Ulnar Rasp, Pilot, Right	

 Table 1 Class II Accessories to the Implant provided non-sterile and are reusable and end-user sterilized



00-8401-034-01	Ulnar Rasp, Size 4/5, Left	
00-8401-034-02	Ulnar Rasp, Size 4/5, Right	
00-8401-036-01	Ulnar Rasp, Size 6, Left	
00-8401-036-02	Ulnar Rasp, Size 6, Right	
00-8401-039-00	Ulnar Clearance Template	
00-8401-044-10	Humeral Provisional, Size 4, 100 mm	
00-8401-044-15	Humeral Provisional, Size 4, 150 mm	
00-8401-045-10	Humeral Provisional, Size 5, 100 mm	
00-8401-045-15	Humeral Provisional, Size 5, 150 mm	
00-8401-046-10	Humeral Provisional, Size 6, 100 mm	
00-8401-046-15	Humeral Provisional, Size 6, 150 mm	
00-8401-055-15	Humeral Provisional, Long Flange, Size 5, 150mm	
00-8401-055-20	Humeral Provisional, Long Flange, Size 5, 200mm	
00-8401-056-15	Humeral Provisional, Long Flange, Size 6, 150mm	
00-8401-056-20	Humeral Provisional, Long Flange, Size 6, 200mm	
00-8401-058-04	Humeral Stem Inserter, Size 4	
00-8401-058-05	Humeral Stem Inserter, Size 5/6	
00-8401-059-04	Humeral Extractor Plate, Size 4	
00-8401-059-05	Humeral Extractor Plate, Size 5/6	
00-8401-060-00	Humeral Awl Reamer	
00-8401-061-00	1/E Alignment Rod	
00-8401-064-00	Humeral Rasp, Pilot	
00-8401-064-10	Humeral Rasp, Size 4, 100 mm	
00-8401-064-15	Humeral Rasp, Size 4, 150 mm	
00-8401-065-10	Humeral Rasp, Size 5, 100 mm	
00-8401-065-15	Humeral Rasp, Size 5, 150 mm	
00-8401-066-10	Humeral Rasp, Size 6, 100 mm	
00-8401-066-15	Humeral Rasp, Size 6, 150 mm	
00-8401-070-45	Flexible Solid Reamer, 4.5 mm	
00-8401-070-05	Flexible Solid Reamer, 5 mm	

 Table 1 (continued) Class II Accessories to the Implant provided non-sterile and are reusable and end-user sterilized



00-8401-072-01	Flexible Cannulated Reamer, 5.5 mm	
00-8401-072-02	Flexible Cannulated Reamer, 6 mm	
00-8401-072-03	Flexible Cannulated Reamer, 6.5 mm	
00-8401-072-04	Flexible Cannulated Reamer, 7 mm	
00-8401-072-05	Flexible Cannulated Reamer, 8 mm	
00-8401-072-06	Flexible Cannulated Reamer, 9 mm	
00-8401-072-07	Flexible Cannulated Reamer, 10 mm	
00-8401-072-08	Flexible Cannulated Reamer, 11 mm	
00-8401-072-09	Flexible Cannulated Reamer, 12 mm	
00-8401-074-00	Trephine, Size 4	
00-8401-075-00	Trephine, Size 5/6	
00-8401-078-00	Humeral Bearing Driver	
00-8401-079-00	Humeral Bearing Driver Pin	
00-8401-081-00	Ulnar Bearing Assembly Tool	
00-8401-082-00	Humeral Bcaring Placement Tool	
00-8401-084-00	Humeral Screw Holder, Size 4	
00-8401-085-00	Humeral Screw Holder, Size 5/6	
00-8401-089-00	Humeral Screw Removal Driver	
00-8401-091-00	Humeral Bearing Extractor	
00-8401-092-00	Articulation Extractor	
00-8401-093-00	Humeral Extractor Screw	
00-8401-097-00	Humeral Screw Removal Driver Bit	
00-8401-004-00	Humeral Cut Guide, Size 4	
00-8401-005-00	Humeral Cut Guide, Size 5/6	

 Table 1 (continued) Class II Accessories to the Implant provided non-sterile and are reusable and end-user sterilized



Of the Class II Accessories to the Implant listed in Section 11a. (Device Description) of the premarket notification submission K123862, the following are provide sterile, are intended to be used during only one surgery (single use) and are not intended to be end-user sterilized (Table 2, below).

Catalog Number	Description	Sterilization Method	Radiation Dose	Validation Method	Sterility Assurance Level (SAL)
00-8401-088-00	Humeral Screw Installation Driver Bit	Gamma Irradiation	25.0-37.0 kGy	AAMI TIR33, VD _{MAX} ²⁰	10 ⁻⁶ or better
00-8401-080-00	Torque-Limiting Humeral Screw Driver	Gamma Irradiation	25.0-40.0 kGy	AAMI 11137-2, VD _{MAX} ²⁵	10 ⁻⁶ or better

Table 2 Class II Accessories to the Implant provided sterile

The Class II Accessories to the Implant (listed in Table 2, above) that are provided sterile are packaged as follows:





Torque-Limiting Humeral Screw Driver

The above listed materials are identical to those currently in use for other similar sterile instrument packages.



FDA Question 4

			Yes	N/A	No			
28.		he device, and/or accessory, and/or a component is provided						
		rile:						
		ect "N/A" if no part of the device, accessories, or						
	<u>co</u> 1	mponents is provided sterile, otherwise complete a-f below.						
	a.	Sterilization method is stated for each component	X	1				
		(including parameters such as dry time for steam						
		sterilization, radiation dose, etc.)		100.00	-0.107			
	b.	A description of method to validate the sterilization cycle	X					
		(e.g., half-cycle method and full citation of FDA-			3			
		recognized standard, including date) is provided						
	c.				x			
		For devices sterilized using chemical sterilants such as						
		ethylene oxide (EO) and hydrogen peroxide, submission						
		states maximum levels of sterilant residuals remaining on						
1	e.	the device and sterilant residual limits.						
		Select "N/A" if not sterilized using chemical sterilants.						
	d.	Submission includes description of packaging and	X					
		packaging contents (e.g., if multiple devices are included						
		within the same package)						
	e.	Sterility Assurance Level (SAL) stated	Х					
	f.				X			
		If device is blood-contacting, a permanent implant, or						
		contacts cerebrospinal fluid, or device is labeled "non-						
		pyrogenic," submission contains a description of the						
		endotoxin method used to make a determination (e.g.,						
		LAL), endotoxin release specification (e.g., 20 EU/device),						
	0	and a rationale for the specification.						
		Select " N/A " if device is not blood-contacting, not a						
		permanent implant, does not contact cerebrospinal fluid,						
		and is not labeled "non-pyrogenic." Select "N/A" if a						
		rationale for omission is provided.						
	Co	Comments: (c) The sponsor has provided a promissory note for EO residuals.						
		Please provide EO residuals remaining on the device following sterilization						
		validation.						
	(f) The sponsor does not provide a rationale for omission of Pyrogen Testin							
		· · · · · · · · · · · · · · · · · · ·						



Zimmer Response to FDA Question 4

- (c) As discussed during the 15 January 2013 phone conversation between Joanna Surma (Zimmer, Inc.) and Peter Allen (FDA), the EO residual data will be provided to Peter Allen via email following each of the following testing phases:
 - EO Residuals after 1 Sterilization Cycle: expected by mid-February
 - EO Residuals after 2 Sterilization Cycles: expected by early March
 - EO Residuals after 3 Sterilization Cycles: expected by early- to mid-March
- (f) Page 87 of the premarket notification submission K123862 stated the following rationale for omission of Pyrogen Testing (additionally, Zimmer makes no "pyrogen-free" or "non-pyrogenic" label claims):

"This device is not labeled as nonpyrogenic. Per USP XXIII (161), requirements for specified endotoxin levels do not apply to orthopaedic implants."



FDA Question 5

		Yes	N/A	No
32.	Submission identifies contact classification (e.g., surface- contacting, less than 24 hour duration)			X
	Comments: Sponsor does not state the contact classification instruments.	ion for	the	

Zimmer Response to FDA Question 5

Of the Class II Accessories to the Implant listed in Section 11a. (Device Description) of the premarket notification submission K123862, all are *external communicating devices* contacting *tissue/bone* with *limited exposure* (less than 24 hour duration), according to definitions presented in Section 5 of ISO 10993-1.

Records processed under F/DIA request #2016-4653; Released by CDRH on 07/16/2017

2/27/13 email



ł

2

Traditional 510(k) Premarket Notification

Predicate Device(s)

The predicate and design basis for the Zimmer Nexel Total Elbow (the proposed device) is the Coonrad/Morrey Total Elbow, manufactured by Zimmer. Clearance letters (K001989 and K053189 cleared 25 July 2000 and 9 December 2005, respectively) for the predicate device are included at the end of this section.

Similarities and differences between the proposed implants and the predicate implants are tabulated below.

Substantial Equivalence Comparison

Indications for Use

	Ēroposed Device: The Zimiñer Nexel Total Elbow	Předičate Device; The Cooprad/Morrey Total Fibow: (K001918, K053189)
Intended Use	To be used in conjunction with bone cement to replace the articulating surfaces of the humeroulnar joint.	To be used in conjunction with bone cement to replace the articulating surfaces of the humeroulnar joint.
Indications for Use	 Elbow joint destruction which significantly compromises the activities of daily living Post-traumatic lesions or bone loss contributing to elbow instability Ankylosed joints, especially in cases of bilateral ankylosis from causes other than active sepsis Advanced rheumatoid, post-traumatic, or degenerative arthritis with incapacitating pain Instability or loss of motion when the degree of joint or soft tissue damage precludes reliable osteosynthesis Acute comminuted articular fracture of the elbow joint surfaces that precludes less radical procedures, including 13-C3 fractures of the distal humerus Revision arthroplasty CAUTION: This device is intended for cemented use only. 	 Elbow joint destruction which significantly compromises the activities of daily living Post-traumatic lesions or bone loss contributing to elbow instability Ankylosed joints, especially in cases of bilateral ankylosis from causes other than sepsis Advanced rheumatoid or degenerative arthritis with incapacitating pain Instability or loss of motion when the degree of joint damage precludes less radical procedures Revision arthroplasty
Product Code	JDC	JDC




Technology

5

	Feature	Proposed Device: The Zimmer Nevel Total Fibow	Predicate Device: The Coonrad/Morrey Total Elbow (K001918, K053189)
	Material	Wrought Tivanium Alloy (Ti-6Al-4V) with Titanium Plasma Spray	Wrought Tivanium Alloy (Ti-6Al-4V) with a porous coating of Titanium Beads
	Fixation Method	Intended for cemented use only	Intended for cemented use only
Humeral Component – General	Comparable Sizes	Size 4, Standard Flange, 100mm $\rightarrow \rightarrow \rightarrow$ Size 4, Standard Flange, 150mm $\rightarrow \rightarrow \rightarrow$ Size 5, Standard Flange, 100mm $\rightarrow \rightarrow \rightarrow$ Size 5, Standard Flange, 150mm $\rightarrow \rightarrow \rightarrow$ Size 5, Long Flange, 150mm $\rightarrow \rightarrow \rightarrow \rightarrow$ Size 6, Standard Flange, 100mm $\rightarrow \rightarrow \rightarrow \rightarrow$ Size 6, Standard Flange, 100mm $\rightarrow \rightarrow \rightarrow \rightarrow$ Size 6, Standard Flange, 150mm $\rightarrow \rightarrow \rightarrow \rightarrow$ Size 6, Long Flange, 150mm $\rightarrow \rightarrow \rightarrow \rightarrow$ Size 6, Long Flange, 200mm $\rightarrow \rightarrow \rightarrow \rightarrow \rightarrow$	 → Size Extra Small, Standard Flange, 4in → Size Extra Small, Standard Flange, 6in → Size Small, Standard Flange, 4in → Size Small, Standard Flange, 6in → Size Small, Long Flange, 6in → Size Regular, Standard Flange, 4in → Size Regular, Standard Flange, 6in → Size Regular, Long Flange, 6in → Size Regular, Long Flange, 6in
	Sterility	Provided Sterile via Gamma Irradiation	Provided Sterile via Gamma Irradiation
	Design Features	 Humeral stem designed to be implanted with bone cement into the patient's humeral medullary canal Anterior flange designed to accept a bone graft and limit torsional and posterior migration Humeral yoke with rounded corners to avoid the creation of stress risers within the medial and lateral humeral supracondylar columns Plasma spray region to enhance fixation to bone cement within the medullary canal, and to improve fatigue strength Features to accept Humeral Screws, Axle Pin and polymer Bearings A and B 	 Humeral stem designed to be implanted with bone cement into the patient's humeral medullary canal Anterior flange designed to accept a bone graft and limit torsional and posterior migration Humeral yoke with squared corners Titanium Bead coating region to enhance fixation to bone cement within the medullary canal Features to accept Axle Pin and polymer Bearings

.








	Feature	Proposed Device: The Zimmer Nevel Total Libow	Predicate Device: The Coonrad/Morrey Total Elbow (K001948, K053189)
	Material	Wrought Tivanium Alloy (Ti-6Al-4V) with Titanium Plasma Spray	Wrought Tivanium Alloy (Ti-6Al-4V) with Titanium Plasma Spray
	Fixation Method	Intended for cemented use only	Intended for cemented use only
Ulnar Component – General	Comparable Sizes	Left, Size 4, 115mm $\rightarrow \rightarrow $	 → Left, Size Extra Small, Standard → Left, Size Extra Small, Long → Left, Size Small, Standard → Left, Size Small, Long → Left, Size Regular, Standard → Left, Size Regular, Long → Right, Size Extra Small, Standard → Right, Size Extra Small, Long → Right, Size Small, Standard → Right, Size Small, Long → Right, Size Small, Long → Right, Size Regular, Standard
L Co	Sterillty	Provided Sterile via Gamma Irradiation	Provided Sterile via Gamma Irradiation
Ulnar	Design Features	 Left and right side-specific Ulnar stem designed to be implanted with bone cement into the patient's ulnar medullary canal Ulnar eye that is both highly polished and nitrogen-enriched to limit wear of the apposing polymer bearings Plasma spray region to enhance fixation to bone cement within the medullary canal 	 Left and right side-specific Ulnar stem designed to be implanted with bone cement into the patient's ulnar medullary canal Ulnar eye that is highly polished to limit wear of the apposing polymer bearings Plasma spray region to enhance fixation to bone cement within the medullary canal








ŝ

Traditional 510(k) Premarket Notification

.

	Feature	Proposed Device: The Zimmer Nevel Total Elbow	Predicate Device: The Coonrad/Morrey Total Elbow (K001918, K053189)
	Material	Vitamin E HXPE (Vitamin E Highly Cross-Linked Polyethylene)	UHMWPE (Ultra-High Molecular-Weight Polyethylene)
	Fixation Method	Press-Fit to Humeral Component	Sliding-Fit to Humeral Component Press-Fit to Ulnar Component
General	Sterility	Provided Sterile via Ethylene Oxide	Provided Sterile via Gamma Irradiation
Polyethylene Components - G	Design Features	 Bearings A and B: Designed to broadly distribute joint reaction forces Designed to articulate against the Eye of the Ulnar Component Features to allow for press-fit mate to Humeral Component Bearings B only: Features to allow for press-fit mate to Axle Pin Designed to articulate against the inner-diameter and the medial and lateral surfaces of the Eye of the Ulnar Component Bearing A only: Designed to articulate against the outer-diameter of the Eye of the Ulnar Component 	 Humeral Bushings: Designed to articulate with the medial and lateral sides of the Eye of the Ulnar Component Features to allow for sliding-fit mate to the Humeral Component and sliding-fit mate to the Axle Pin (also known as Hinge Pin, also known as Snap Pin) Ulnar Bushing: Features to allow for a press-fit mate with the Eye of the Ulnar Component Designed to articulate against the Humeral Bushings







	Feature	Proposed Device: The Zimmer Nevel Total Elbow	Prédicate Device: The Coonrad/Morrey Total Elbow (K001918, K053189)
	Material	Zimaloy (Cobalt-Chromium-Molybdenum Alloy)	Zimaloy (Cobalt-Chromium-Molybdenum Alloy)
Axle Pin	Fixation Method	 Press-Fit to Bearings B Compression fit between the taper of the Humeral Screws and the yoke of the Humeral Component 	 Sliding-Fit to Bearings Snap-Fit between inner and outer Axle Pin components
	Sizes	Size 4 (b) Long Diameter • Solid Material Throughout Size 5/6 (b) Long Diameter • Solid Material Throughout	Size Extra Small (b) (4) "Outer Diameter "Wall Thickness Size Small/Regular (b)(4) Long Outer Diameter Wall Thickness
	Sterility	Provided Sterile via Ethylene Oxide	Provided Sterile via Gamma Irradiation
	; Design Features	 Designed to remain in a fixed position (not translate nor rotate) throughout the entire humeroulnar joint range of motion Machined complete as a single, solid component 	 Designed to freely rotate and slide Machined as 2 (an inner and an outer) hollow components that snap together intraoperatively
	Material	Zimaloy (Cobalt-Chromium-Molybdenum Alloy)	
ľ	Fixation Method	Threads engage with the Spiralock® threads of the Humeral Component	
	Sizes	Single size mates with all humeral component sizes and both axle pin sizes	
crews	Sterility	Provided Sterile via Gamma Irradiation	
Humeral Screws	Design Features	 Designed to remain in a fixed position throughout the entire humeroulnar joint range of motion (not directly loaded by any compressive joint reaction force during flexion or extension of the humeroulnar joint) Headless design intended to ensure the constant application of a clamp-load to the axle-pin and to avoid soft tissue impingement or irritation associated with a protruding screw head Non-threaded taper to (along with the Humeral Component) hold Axle Pin in compression Pilot feature to avoid cross-threading 	Not Applicable: No comparable part in the predicate system



b,

		Constraint Characteristics	Próposed Device: The Zimmer Nevel Total Elbow	Predicate Device: The Coonrad/Morrey Total Elbow (K001918, K053189)
mponent)		Varus/Valgus Angular Laxity	Neutral to Maximum Varus: (b) Neutral to Maximum Valgus: (b) Total Included Angle: (b)(Neutral to Maximum Varus: (b)(Neutral to Maximum Valgus: 4755 Total Included Angle: (b)
to Humeral Co	Internal/External Rotational Laxity		Neutral to Maximum Internal Rotation: (b)(Neutral to Maximum External Rotation: (b)(Total Included Angle: (b)	Neutral to Maximum Internal Rotation: (b) Neutral to Maximum External Rotation: (b) Total Included Angle: (b)
tion respect	Pistoning	Proximal/Distal Translational Laxity	Neutral to Maximum Proximal (b)(4) Neutral to Maximum Distal: Total Proximal/Distal Laxity:	Neutral to Maximum Proximal: (b)(4) Neutral to Maximum Distal: Total Proximal/Distal Laxity:
ge of Mo ent with	Pisto	Anterior/Posterior Translation Laxity	Neutral to Maximum Anterior: Neutral to Maximum Posterior: Total Anterior/Posterior Laxity	Neutral to Maximum Anterior: Neutral to Maximum Posterior: Total Anterior/Posterior Laxity:
ned Ran t movem		Medial/Lateral Translation Laxity	Neutral to Maximum Medial: Neutral to Maximum Lateral: Total Medial/Lateral Laxity:	Neutral to Maximum Medial: Neutral to Maximum Lateral: Total Medial/Lateral Laxity:
Constrained Range of Motion (defined in terms of Ulnar Component movement with respect to Humeral Component)		Flexion/Extension Range defined in terms of degrees flexion)	All Humeral and Ulnar combinations: (b)(4)	Humeral size XS with Ulnar size XS: (b)(4) Humeral size XS with Ulnar size S: (c)(4) Humeral size XS with Ulnar size REG: (b)(4) Humeral size S with Ulnar size XS: (b)(4) Humeral size S with Ulnar size S: (b)(4) Humeral size REG with Ulnar size XS: (b)(4) Humeral size REG with Ulnar size S: (b)(4) Humeral size REG with Ulnar size S: (b)(4) Humeral size REG with Ulnar size REG: (b)(4) Humeral size REG with Ulnar size REG: (b)(4) Humeral size REG with Ulnar size REG: (b)(4)



Why Any Differences Do Not Render the Device NSE

The major technological differences between the proposed device (Zimmer Nexel Total Elbow) and the predicate device (the Zimmer Coonrad/Morrey Total Elbow) are listed below along with the performance tests conducted to demonstrate substantial equivalence (the listed performance tests and their conclusions are described in more detail in *Section 18a Performance Testing – Bench* of this premarket notification submission):

Device Component	Major Difference Between Proposed and Predicate Devices	Performance Test Conducted to Demonstrate Substantial Equivalence
Humeral Component	 Distal design features for implant assembly Plasma Spray 	 Durability Testing Modular Connection Fatigue Testing Fatigue Testing
Ulnar Component	- Nitrogen enrichment	 Wear Testing Fatigue Testing
Polyethylene Components	 Polyethylene material Surfaces against which components are intended to articulate 	Wear TestingDurability Testing
Axle Pin	- Single Solid vs. 2-Part Hollow construction	 Wear Testing Durability Testing Modular Connection Fatigue Testing
Humeral Screws	- Existence of screws	 Wear Testing Durability Testing Modular Connection Fatigue Testing

Current questions of safety and effectiveness (clinically observed failure modes) associated with elbow prostheses are humeral and ulnar component loosening, humeral and ulnar component fracture, polyethylene wear, system (implant assembly) durability and system (implant assembly) disassociation (also known as modular connection fatigue). To demonstrate that the proposed device is at least as safe and effective as the predicate device despite the compounding effects of these technological differences and previously outlined dimensional differences, fatigue (to assess resistance to fracture), wear, durability and modular connection fatigue performance tests were conducted (testing conclusions are described in more detail in *Section 18a Performance*



Testing – Bench of this premarket notification submission). Stem loosening was not the focus of a performance test for the New Zimmer Total Elbow, because the variables (cementing technique, stem geometry, stem material and stem surface finish) associated with stem loosening are consistent from the predicate to the proposed device; and, therefore the proposed device (the New Zimmer Total Elbow) is expected to be as safe and effective as the predicate device (the Coonrad/Morrey Total Elbow) in terms of ability to resist stem loosening.

As reflected in the dFMEA documentation associated with the proposed device, the Zimmer Nexel Total Elbow does not raise any new/different types of questions of safety or effectiveness.

The proposed and predicate devices have the exact same intended use. The proposed and predicate devices share similar indications for use; where, the only changes in wording of the indications for use are that the proposed device has indications for use that are intended to add specificity to the indications for use of the proposed device. This added specificity does not change the intended meaning of the indications for use of the proposed device as compared to those of the predicate.

Because the aforementioned differences between the proposed and predicate devices do not constitute a new intended use, affect safety or effectiveness or raise different questions of safety and effectiveness, these differences *do not* render the proposed device NSE (not substantially equivalent).



.

Traditional 510(k) Premarket Notification

Performance Specifications

	Test Resúlts			
Test	Proposed Device: The Zimmer Nexel Total Elbow	Predicate Device: The Coonrad/Morrey Total Elbow (K001918, K053189)		
Fatigue Testing: Humeral Stem (Plasma Spray Region)	Five Zimmer Nexel Total Elbow humeral stems achieved (b)(4) load cycles at (b)(4) without fracture.	The predicate device humeral component was fatigue tested at various load levels, and the fatigue strength at $(6)(4)$ cycles was determined to be $755/4$		
Fatigue Testing: Uluar Stem (Mid-Stem Region)Five Zimmer Nexel Total Elbow ulna stems achieved (b)(4) total Elbow ulna 		The predicate device ulnar component was fatigue tested at various load levels, and the fatigue strength at (h)(4) toad cycles in the mid-stem region was determined to be (h)(4)		
Fatigue Testing: Ulnar Stem (Plasma Spray Region)	Five Zimmer Nexel Total Elbow ulnar stems achieved (b)(4). I load cycles in the plasma spray region at a load of (b)(4) without fracture.	The predicate device ulnar component was fatigue tested at various load levels, and the fatigue strength at (b)(4). I load cycles in the plasma spray region was determined to be (b)(
Wear Testing	The Zimmer Nexel Total Elbow mean gravimetric wear rate $(b)(4)$ at $(b)(4)$ load cycles) was 91% less than that of the predicate.	The predicate device mean gravimetric wear rate was $\frac{(b)(4)}{(b)(4)}$ load cycles.		
Durability Testing	The Zimmer Nexel Total Elbow device achieved (h)(4) run-outs at an equivalent of (b)(weight-in-hand.	The run out load for the predicate device was found to be equivalent to (b) bf weight-in-hand.		
Modular Connection Fatigue Testing	The median fatigue strength of the non- articulating, mechanically locked, modular Zimmer Nexel Total Elbow implant components was no less than (b)(4).	Not Applicable: The Zimmer Coonrad/Morrey Total Elbow does not have non-articulating, mechanically-locked, modular implant components comparable to those of the proposed device.		

2

.

Allen, Peter

From:	Jason Heckaman <jason.heckaman@zimmer.com></jason.heckaman@zimmer.com>
Sent:	Thursday, Eebruary 14, 2013 12:55 PM
To:	Allen, Peter
Cc:	Joanna L. Surma
Subject:	RE: K123862
Attachments:	#2013-02-14-Nexel-Articulation Kit EO Residual-Results:pdf-

Hi Pete,

Please find enclosed EO residual data at 16 hours, 2-day, and 3-day time points following both one and two sterilization cycles. All results were within the specifications.

If there is anything else I can help you with, just let me know.

Kind Regards, Jason

.

From: Jason Heckaman Sent: Tuesday, February 12, 2013 10:58 AM To: 'Allen, Peter' Cc: Joanna L. Surma Subject: RE: K123862

Hi Pete,

By the end of business Friday we expect to be able to provide you with EO residual data following both one and two sterilization cycles.

Kind Regards, Jason

From: Jason Heckaman Sent: Monday, February 11, 2013 5:10 PM To: 'Allen, Peter'; Joanna L. Surma Subject: RE: K123862

Hi Pete,

I will come back to you with an update by the end of the day tomorrow.

Kind Regards, Jason

From: Allen, Peter [<u>mailto:Peter,Allen@fda.hhs.gov]</u> Sent: Monday, February 11; 2013 2:21 PM To: Joanna L. Surma Cc: Jason Heckaman Subject: RE: K123862

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

1

Joanna/Jason,

I see that Joanna is out of the office this week, so Jason I am hoping you will be able to respond to this email ASAP. This is in regards to 510(k) submission K123862 – Zimmer Nexel Total Elbow.

Joanna had discussed the possibility of getting Ethylene Oxide residual testing results to me by mid-February. I need to know if these results will be available in the next week. Otherwise I will place this submission on hold until the results are available. Please let me know if you will be able to accommodate my request.

Thanks,

Pete

Peter G. Allen, MS Biomedical Engineer Joint Fixation One Branch (JFOB) FDA/CDRH/ODE/DOD 10903 New Hampshire Avenue Silver Spring, MD 20993-0002 Ph: (301) 796-6402 Fax: (301) 847-8119 Email: <u>peter.allen@fda.hhs.gov</u>

This e-mail message is intended for the exclusive use of the recipient(s) named above. It may contain information that is protected, privileged, or confidential, and it should not be disseminated, distributed, or copied to persons not authorized to receive such information. If you are not the intended recipient, any dissemination, distribution or copying is strictly prohibited. If you think you have received this e-mail message in error, please e-mail the sender immediately at peter.allen@fda.hhs.goy

From: Joanna L. Surma [<u>mailto:Joanna.Surma@zimmer.com</u>] Sent: Friday, February 08, 2013 9:19 AM To: Allen, Peter Cc: Jason Heckaman Subject: K123862

Dear Mr. Allen,

I will be out of the office next week (10-15 Feb). Could you please copy my manager, Jason Heckaman, on any communications to me regarding K123862 (Zimmer Nexel Total Elbow) during that time?

Jason's email is Jason. Heckaman@Zimmer.com, and he can be reached by phone at 574.371.8675.

Thank you, Joanna

Joanna L. Surma - MS, RAC Associate Project Manager Regulatory Affairs Zimmer, Inc. (o) 574.371.1642

2

EO Residual Results

Zimmer Nexel™ Total Elbow Articulation Kit

Cycle + Timepoint	Analyte	Specification	Result
1X @ 16 hours	Ethylene	≤4mg/1 st 24 hours	(b)(4)
	Oxide	≤60mg/1 st 30 days	
	(EO)	≤2.5g lifetime	
		≤0.1 mg/day ADD*	
	Ethylene	≤9mg/1 st 24 hours	
	Chlorohydrin	≤60mg/1 st 30 days	
	(ECH)	≤10g lifetime	
		≤0.4 mg/day ADD	
	TCL** EO	≤10µg/cm²	
	TCL ECH	≤5mg/cm ²	
1X @ Day 2	Ethylene	≤4mg/1 st 24 hours	
	Oxide	≤60mg/1 st 30 days	
	(EO)	≤2.5g lifetime	
		≤0.1 mg/day ADD	
	Ethylene	≤9mg/1 st 24 hours	
	Chlorohydrin	≤60mg/1 st 30 days	
	(ECH)	≤10g lifetime	
		≤0.4 mg/day ADD	
	TCL EO	≤10µg/cm²	
	TCL ECH	≤5mg/cm ²	
1X @ Day 3	Ethylene	≤4mg/1 st 24 hours	
	Oxide	≤60mg/1 st 30 days	
	(EO)	≤2.5g lifetime	
		≤0.1 mg/day ADD	
	Ethylene	≤9mg/1 st 24 hours	
	Chlorohydrin	≤60mg/1 st 30 days	
	(ECH)	≤10g lifetime	
		≤0.4 mg/day ADD	
	TCL EO	≤10µg/cm²	
	TCL ECH	≤5mg/cm ²	
*Average Daily Dose		Confact Limit	

*<u>Average Daily Dose</u> **<u>Tolerable Contact Limit</u>

2

1

Cycle + Timepoint	Analyte	Specification	Result
2X @ 16 hours	Ethylene	≤4mg/1 st 24 hours	(b)(4)
100-5	Oxide	≤60mg/1 st 30 days	
	(EO)	≤2.5g lifetime	
		≤0.1 mg/day ADD*	
	Ethylene	≤9mg/1 st 24 hours	
	Chlorohydrin	≤60mg/1 st 30 days	
	(ÉCH)	≤10g lifetime	
		≤0.4 mg/day ADD	
	TCL** EO	≤10µg/cm²	
	TCL ECH	≤5mg/cm ²	
2X @ Day 2	Ethylene	≤4mg/1 st 24 hours	
	Oxide	≤60mg/1 st 30 days	
	(EO)	≤2.5g lifetime	
		≤0.1 mg/day ADD	
	Ethylene	≤9mg/1 st 24 hours	
	Chlorohydrin	≤60mg/1 st 30 days	
	(ECH)	≤10g lifetime	
		≤0.4 mg/day ADD	
	TCL EO	≤10µg/cm²	
	TCL ECH	≤5mg/cm ²	
2X @ Day 3	Ethylene	≤4mg/1 st 24 hours	
	Oxide	≤60mg/1 st 30 days	
	(EO)	≤2.5g lifetime	
		≤0.1 mg/day ADD	
	Ethylene	≤9mg/1 st 24 hours	
	Chlorohydrin	≤60mg/1 st 30 days	
	(ECH)	≤10g lifetime	
		≤0.4 mg/day ADD	
	TCL EO	≤10µg/cm ²	
	TCL ECH	≤5mg/cm ²	
*Average Daily Dose	** <u>T</u> olerable	<u>C</u> ontact <u>L</u> imit	

Allen, Peter

From:	Joanna L. Surma <joanna.surma@zimmer.com></joanna.surma@zimmer.com>
Sent:	Wednesday, aniuary 16, 2013 11:05 AM
To:	Allen, Peter
Subject:	RE: K123862 Refuse to Accept Notification - Zimmer Response
₽Attachments: ?:	C=20120101-2013-01-17 002.01 RTA Response.pdf

Dear Mr. Allen,

Thank you for taking the time to call me yesterday to discuss the EO residual data availability timeline.

Within the attached RTA response, you'll find the following summary of that conversation:

As discussed during the 15 January 2013 phone conversation between Joanna Surma (Zimmer, Inc.) and Peter Allen (FDA), the EO residual data will be provided to Peter Allen via email following each of the following testing phases:

- EO Residuals after 1 Sterilization Cycle: expected by mid-February
- EO Residuals after 2 Sterilization Cycles: expected by early March
- EO Residuals after 3 Sterilization Cycles: expected by early- to mid-March

If you require any additional information or have any questions, please contact me by telephone (574.371.1642), email (Joanna.Surma@Zimmer) or fax (574.372.4605).

Sincerely, Joanna

Joanna L. Surma - MS, RAC Associate Project Manager Regulatory Affairs Zimmer, Inc. (o) 574.371.1642 (c) 574.527.2361

From: Joanna L. Surma Sent: Monday, January 14, 2013 8:40 PM To: Allen, Peter Cc: 'Hanley, Casey' Subject: RE: K123862 - Refuse to Accept Notification

Dear Mr. Allen,

Thank you and Dr. Hanley for your feedback.

We have completed our response to your initial review of our submission for completeness, and would like to discuss the acceptability of a proposed timeline of availability for EO residuals (Acceptance Checklist Item 28c) via a 10-15minute phone conference before submitting the formal response. Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017

Please let me know you're availability Tuesday through Thursday this week for that phone conference, as we intend to send the formal response by Thursday 17 January 2013 to ensure its delivery to the FDA by 21 January 2013.

Sincerely, Joanna

Joanna L. Surma - MS, RAC Associate Project Manager Regulatory Affairs Zimmer, Inc. (o) 574.371.1642 (c) 574.527.2361

From: Hanley, Casey [<u>mailto:Casey.Hanley@fda.hhs.gov</u>] Sent: Friday, December 28, 2012 3:49 PM To: Joanna L. Surma Cc: Allen, Peter Subject: K123862 - Refuse to Accept Notification

Dear Ms. Surma,

We have completed the administrative acceptance review of your premarket notification (510(k)) submission K123862) Our review indicates that your 510(k) submission does not meet the criteria established for administrative completeness.

Please refer to the attached checklist and e-mail your response referencing the 510(k) number K123862 by January 21, 2013. Your response should address all of the elements identified as missing or inconsistent in the attached checklist. Failure to provide the requested information by January 21, 2013 may result in placing the file on hold.

Upon receipt of the requested information, FDA may have additional requests for information.

Please be advised that once our RTA policy is implemented, if your 510(k) submission does not meet the criteria established for administrative completeness, your file will not be accepted and we will not begin our substantive review until you submit the missing elements and your submission is accepted. Refer to the draft <u>guidance</u> document for information.

Should you have questions about this email, you may either contact me or Peter Allen, the lead reviewer assigned to your 510(k) submission.

Sincerely, Casey

Casey L. Hanley, Ph.D. U.S. Food and Drug Administration CDRH/ODE/DOD/JFOB 10903 New Hampshire Ave WO66 Room 1567 Silver Spring, MD 20993 Tel: (301)796 6948 <u>casey.hanley@fda.hhs.gov</u>

This e mail message is intended for the exclusive use of the recipient(s) named above. It may contain information that is protected, privileged, or confidential and it should not be disseminated distributed or copied to persons not authorized to receive such information. If you are not the intended recipient, any dissemination distribution or copying is strictly prohibited if you think you have received this e mail message in error please e-mail the sender immediately.

2

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017

This communication is consistent with 21 CFR 10.85 (k) and constitutes an informal communication that represents my best judgment at this time but does not constitute an advisory opinion, does not necessarily represent the formal position of FDA, and does not bind or otherwise obligate or commit the agency to the views expressed.

Table of Contents

Submission	2
Cover Letter	2
Zimmer Responses	3
FDA Question 1	3
Zimmer Response to FDA Question 1	4
FDA Question 2	5
Zimmer Response to FDA Question 2	6
FDA Question 3	8
Zimmer Response to FDA Question 3	9
FDA Question 4	13
Zimmer Response to FDA Question 4	14
FDA Question 5	15
Zimmer Response to FDA Question 5	15
Attachment 01 - 510(k) Summary	16
Attachment 02 - Substantial Equivalence Discussion	19
Attachment 03 - Cleaning and ReSterilization Instructions	31

L



P.O. Box 708 Warsaw, IN 46581-0708 (574) 267-6131

17 January 2013

U.S. Food and Drug Administration Center for Devices and Radiological Health Document Mail Center – W066-G609 10903 New Hampshire Avenue Silver Spring, MD 20993-0002

Dear Sir or Madam:

Subject: Additional Information – Zimmer Nexel Total Elbow – 510(k) (K123862)

Enclosed are Zimmer's responses to FDA's questions regarding the above-referenced submission.

If you require any additional information or have any questions, please contact me by telephone at (574) 371-1642, by e-mail at Joanna.Surma@Zimmer.com or by fax at (574) 372-4605.

Sincerely,

1 Jan 2013 O'L Sum

Joanna L. Surma Associate Project Manager, Regulatory Affairs

Enclosure

ĩ.



FDA Question 1

0			Yes	N/A	No	
4.	Eit	bmission contains 510(k) Summary or 510(k) Statement her a) or b) must be answered "Yes" to be considered nplete. Identify any missing element(s) as Comments.	X			
	a.	Summary contains all elements per 21 CFR 807.92 See also 510(k) Summary Checklist			X	
	b .	Statement contains all elements per 21 CFR 807.93	6			
	Comments: The 510(k) summary is missing a complete Device Description, (lacking explanation of how the device functions, the scientific concepts that form the basis for the device, and the significant physical and performance characteristics of the device, such as device design, material used, and physical properties) and is missing a comparison to predicate devices.					

CONFIDENTIAL AND PROPRIETARY INFORMATION EXEMPT FROM DISCLOSURE UNDER FOI



Zimmer Response to FDA Question 1

Please see Attachment 1 for a replacement 510(k) summary.

Within the device description section of the replacement 510(k) summary, the following verbiage was added to that of the original submission:

How the Device Functions: The Zimmer Nexel Total Elbow is an implant designed to replace the articulating surfaces of and restore motion to the humeroulnar (elbow) joint. The implant is a semiconstrained device assembly and consists of the following components: humeral component, ulnar component, humeral bearing-A, ulnar bearings-B, axle pin and humeral screws.

The Following Scientific Concepts, Design Features and Physical Properties form the Basis for the Zimmer Nexel Total Elbow: The humeral component has a humeral stem designed to be implanted with bone cement into the patient's humeral medullary canal, an anterior flange designed to accept a bone graft and limit torsional and posterior migration, a humeral yoke with rounded corners to avoid the creation of stress risers within the medial and lateral humeral supracondylar columns, and plasma spray region to enhance fixation to bone cement within the medullary canal, and to improve fatigue strength. The ulnar component has an ulnar stem designed to be implanted with bone cement into the patient's ulnar medullary canal, an ulnar eye that is both highly polished and nitrogen-enriched to limit wear of the apposing polymer bearings, and plasma spray region to enhance fixation to bone cement within the medullary canal. Bearings A and B are designed to broadly distribute joint reaction forces.

Materials Used: The humeral and ulnar components are made of a titanium alloy, the bearings A and B are made of Vitamin E highly cross-linked polyethylene, and the humeral screws and axle pin are made of a cobalt-chromium-molybdenum alloy.

Comparison to the Predicate: The proposed device (Zimmer Nexel Total Elbow) and the predicate device (Coonrad/Morrey Total Elbow) have the same intended use and similar indications for use. The proposed device humeral and ulnar components are very similar in terms of materials used and design/dimensions to the predicate device humeral and ulnar components. The primary differences between the proposed and predicate devices are that the proposed device bearing components are made of Vitamin E highly cross-linked polyethylene (the predicate device bearing components are made of ultrahigh molecular-weight polyethylene), and the proposed device has bearings that articulate on both the outer and inner diameters of the ulnar eye (the predicate device uses screws and a 1-part axle pin to fix the bearings in place, while the predicate device uses a 2-part (snap-fit) axle pin.

CONFIDENTIAL AND PROPRIETARY INFORMATION EXEMPT FROM DISCLOSURE UNDER FOI

4



FDA Question 2

		Yes	N/A	No
17.	Submission includes an analysis of why any differences between the subject device and predicate(s) do not render the device NSE (e.g., do not constitute a new intended use, affect safety or effectiveness, or raise different questions of safety and effectiveness) (see section 513(i)(1)(A) of the FD&C Act)			x
	If there is no difference between the subject and predicate(s) with respect to indications for use or technology, this should be explicitly stated, in which case "N/A" should be selected. Select "No" only if the submission does not include an analysis of differences as described above or a statement that there are no differences. Note that due to potential differences in manufacturing that may not be known to the submitter, no identified differences does not necessarily mean that no performance testing is needed.			•
	Comments: The sponsor provides a tabulated comparison of and predicate devices. From this tabulated list, it is clear the differences between the subject and predicate devices with r indications for use and design features. The sponsor does no analysis/discussion of why these differences do not render the	ere are espec of prov	e t to vide ar	

CONFIDENTIAL AND PROPRIETARY INFORMATION EXEMPT FROM DISCLOSURE UNDER FOI



Zimmer Response to FDA Question 2

Please see <u>Attachment 2</u> for a replacement Substantial Equivalence Discussion.

Within the replacement Substantial Equivalence Discussion, the following verbiage was added to address FDA Question 2:

Why Any Differences Do Not Render the Device NSE

The major technological differences between the proposed device (Zimmer Nexel Total Elbow) and the predicate device (the Zimmer Coonrad/Morrey Total Elbow) are listed below along with the performance tests conducted to demonstrate substantial equivalence (the listed performance tests and their conclusions are described in more detail in *Section 18a Performance Testing – Bench* of this premarket notification submission):

Device Component	Major Difference Between Proposed and Predicate Devices	Performance Test Conducted to Demonstrate Substantial Equivalence
Humeral Component	 Distal design features for implant assembly Plasma Spray 	 Durability Testing (ZRR_WA_2542_12) Modular Connection Fatigue Testing (ZRR_WA_2598_12) Fatigue Testing (ZRR_WA_2554_12 & ZRR_WA_2589_12)
Ulnar Component	- Nitrogen enrichment	 Wear Testing (ZRR_WA_2407_11 & ZRR_WA_2552_12) Fatigue Testing (ZRR_WA_2585_12, ZRR_WA_2615_12, ZRR_WA_2554_12 & ZRR_WA_2626_12)
Polyethylene Components	 Polyethylene material Surfaces against which components are intended to articulate 	 Wear Testing (ZRR_WA_2407_11 & ZRR_WA_2552_12) Durability Testing (ZRR_WA_2542_12)
Axle Pin	- Single Solid vs. 2-Part Hollow construction	 Wear Testing (ZRR_WA_2407_11 & ZRR_WA_2552_12) Durability Testing (ZRR_WA_2542_12) Modular Connection Fatigue Testing (ZRR_WA_2598_12)
Humeral Screws	- Existence of screws	 Wear Testing (ZRR_WA_2407_11 & ZRR_WA_2552_12) Durability Testing (ZRR_WA_2542_12) Modular Connection Fatigue Testing (ZRR_WA_2598_12)



Current questions of safety and effectiveness (clinically observed failure modes) associated with elbow prostheses are humeral and ulnar component loosening, humeral and ulnar component fracture, polyethylene wear, system (implant assembly) durability and system (implant assembly) disassociation (also known as modular connection fatigue). To demonstrate that the proposed device is at least as safe and effective as the predicate device despite the compounding effects of these technological differences and previously outlined dimensional differences, fatigue (to assess resistance to fracture), wear, durability and modular connection fatigue performance tests were conducted (testing conclusions are described in more detail in *Section 18a Performance Testing – Bench* of this premarket notification submission). Stem loosening was not the focus of a performance test for the New Zimmer Total Elbow, because the variables (cementing technique, stem geometry, stem material and stem surface finish) associated with stem loosening are consistent from the predicate to the proposed device; and, therefore the proposed device (the New Zimmer Total Elbow) is expected to be as safe and effective as the predicate device (the Coonrad/Morrey Total Elbow) in terms of ability to resist stem loosening.

As reflected in the dFMEA documentation associated with the proposed device, the Zimmer Nexel Total Elbow does not raise any new/different types of questions of safety or effectiveness.

The proposed and predicate devices have the exact same intended use. The proposed and predicate devices share similar indications for use; where, the only changes in wording of the indications for use are that the proposed device has indications for use that are intended to add specificity to the indications for use of the proposed device. This added specificity does not change the intended meaning of the indications for use of the proposed device as compared to those of the predicate.

Because the aforementioned differences between the proposed and predicate devices do not constitute a new intended use, affect safety or effectiveness or raise different questions of safety and effectiveness, these differences *do not* render the proposed device NSE (not substantially equivalent).



FDA Question 3

		· · · · · · · · · · · · · · · · · · ·	Yes	N/A	No
27.	Assessment of the need for sterilization information				
	a.	Identification of device, and/or accessories, and/or components that are provided sterile.	x		
	b.	Identification of device, and/or accessories, and/or components that are end user sterilized			X
	ins pro re- cor	mments: The sponsor has identified a number of patient of truments; however, the sponsor has not discussed whethe ovided sterile or non-sterile and whether they are able to sterilized. Therefore, please provide 1) a discussion as to nponents are end user sterilized and 2) instructions for cl rilization for these components.	er they be clear which	are ned an	



Zimmer Response to FDA Question 3

Of the Class II Accessories to the Implant listed in Section 11a. (Device Description) of the premarket notification submission K123862, the following are provided non-sterile and are reusable and end-user sterilized (Table 1, below). Instructions for cleaning and re-sterilization of these instruments are provided in <u>Attachment 3</u>.

Catalog Number	Description	
00-8401-002-00	T-Handle	
00-8401-008-00	Reciprocating Saw Adapter	
00-8401-008-01	Single Use Connector	
00-8401-009-00	Slide Hammer	
00-8401-018-00	Ulnar Bearing Tamp	
00-8401-012-04	Trephine Stabilizer, Size 4	
00-8401-012-05	Trephine Stabilizer, Size 5	
00-8401-014-07	Ulnar Provisional, Size 4, 75 mm, Left	
00-8401-014-11	Ulnar Provisional, Size 4, 115 mm, Left	
00-8401-015-07	Ulnar Provisional, Size 5, 75 mm, Left	
00-8401-015-11	Ulnar Provisional, Size 5, 115 mm, Left	
00-8401-016-09	Ulnar Provisional, Size 6, 90 mm, Left	
00-8401-016-11	Ulnar Provisional, Size 6, 115 mm, Left	
00-8401-024-07	Ulnar Provisional, Size 4, 75 mm, Right	
00-8401-024-11	Ulnar Provisional, Size 4, 115 mm, Right	
00-8401-025-07	Ulnar Provisional, Size 5, 75 mm, Right	
00-8401-025-11	Ulnar Provisional, Size 5, 115 mm, Right	
00-8401-026-09	Ulnar Provisional, Size 6, 90 mm, Right	
00-8401-026-11	Ulnar Provisional, Size 6, 115 mm, Right	
00-8401-019-00	Articulation Inserter	
00-8401-028-00	Ulnar Stem Inserter	
00-8401-029-00	Implant Extractor Hook	
00-8401-033-01	Ulnar Rasp, Pilot, Left	
00-8401-033-02	Ulnar Rasp, Pilot, Right	

 Table 1 Class II Accessories to the Implant provided non-sterile and are reusable and end-user sterilized

CONFIDENTIAL AND PROPRIETARY INFORMATION EXEMPT FROM DISCLOSURE UNDER FOI

9

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017



K123862 Zimmer Nexel Total Elbow Refuse to Accept – Zimmer Response

00-8401-034-01	Ulnar Rasp, Size 4/5, Left	
00-8401-034-02	Ulnar Rasp, Size 4/5, Right	
00-8401-036-01	Ulnar Rasp, Size 6, Left	
00-8401-036-02	Ulnar Rasp, Size 6, Right	
00-8401-039-00	Ulnar Clearance Template	
00-8401-044-10	Humeral Provisional, Size 4, 100 mm	
00-8401-044-15	Humeral Provisional, Size 4, 150 mm	
00-8401-045-10	Humeral Provisional, Size 5, 100 mm	
00-8401-045-15	Humeral Provisional, Size 5, 150 mm	
00-8401-046-10	Humeral Provisional, Size 6, 100 mm	
00-8401-046-15	Humeral Provisional, Size 6, 150 mm	
00-8401-055-15	Humeral Provisional, Long Flange, Size 5, 150mm	
00-8401-055-20	Humeral Provisional, Long Flange, Size 5, 200mm	
00-8401-056-15	Humeral Provisional, Long Flange, Size 6, 150mm	
00-8401-056-20	Humeral Provisional, Long Flange, Size 6, 200mm	
00-8401-058-04	Humeral Stem Inserter, Size 4	
00-8401-058-05	Humeral Stem Inserter, Size 5/6	
00-8401-059-04	Humeral Extractor Plate, Size 4	
00-8401-059-05	Humeral Extractor Plate, Size 5/6	
00-8401-060-00	Humeral Awl Reamer	
00-8401-061-00	I/E Alignment Rod	
00-8401-064-00	Humeral Rasp, Pilot	
00-8401-064-10	Humeral Rasp, Size 4, 100 mm	
00-8401-064-15	Humeral Rasp, Size 4, 150 mm	
00-8401-065-10	Humeral Rasp, Size 5, 100 mm	
00-8401-065-15	Humeral Rasp, Size 5, 150 mm	
00-8401-066-10	Humeral Rasp, Size 6, 100 mm	
00-8401-066-15	Humeral Rasp, Size 6, 150 mm	
00-8401-070-45	Flexible Solid Reamer, 4.5 mm	
00-8401-070-05	Flexible Solid Reamer, 5 mm	

 Table 1 (continued) Class II Accessories to the Implant provided non-sterile and are reusable and end-user sterilized

CONFIDENTIAL AND PROPRIETARY INFORMATION EXEMPT FROM DISCLOSURE UNDER FOI



00-8401-072-01	Flexible Cannulated Reamer, 5.5 mm
00-8401-072-02	Flexible Cannulated Reamer, 6 mm
00-8401-072-03	Flexible Cannulated Reamer, 6.5 mm
00-8401-072-04	Flexible Cannulated Reamer, 7 mm
00-8401-072-05	Flexible Cannulated Reamer, 8 mm
00-8401-072-06	Flexible Cannulated Reamer, 9 mm
00-8401-072-07	Flexible Cannulated Reamer, 10 mm
00-8401-072-08	Flexible Cannulated Reamer, 11 mm
00-8401-072-09	Flexible Cannulated Reamer, 12 mm
00-8401-074-00	Trephine, Size 4
00-8401-075-00	Trephine, Size 5/6
00-8401-078-00	Humeral Bearing Driver
00-8401-079-00	Humeral Bearing Driver Pin
00-8401-081-00	Ulnar Bearing Assembly Tool
00-8401-082-00	Humeral Bearing Placement Tool
00-8401-084-00	Humeral Screw Holder, Size 4
00-8401-085-00	Humeral Screw Holder, Size 5/6
00-8401-089-00	Humeral Screw Removal Driver
00-8401-091-00	Humeral Bearing Extractor
00-8401-092-00	Articulation Extractor
00-8401-093-00	Humeral Extractor Screw
00-8401-097-00	Humeral Screw Removal Driver Bit
00-8401-004-00	Humeral Cut Guide, Size 4
00-8401-005-00	Humeral Cut Guide, Size 5/6

 Table 1 (continued) Class II Accessories to the Implant provided non-sterile and are reusable and end-user sterilized

CONFIDENTIAL AND PROPRIETARY INFORMATION EXEMPT FROM DISCLOSURE UNDER FOI



Of the Class II Accessories to the Implant listed in Section 11a. (Device Description) of the premarket notification submission K123862, the following are provide sterile, are intended to be used during only one surgery (single use) and are not intended to be end-user sterilized (Table 2, below).

Catalog Number	Description	Sterilization Method	Radiation Dose	Validation Method	Sterility Assurance Level (SAL)
00-8401-088-00	Humeral Screw Installation Driver Bit	Gamma Irradiation	(b)(4)	AAMI TIR33, VD _{MAX} ²⁰	10 ⁻⁶ or better
00-8401-080-00	Torque-Limiting Humeral Screw Driver	Gamma Irradiation		AAMI 11137-2, VD _{MAX} ²⁵	10 ⁻⁶ or better


Table 2 Class II Accessories to the Implant provided sterile

The Class II Accessories to the Implant (listed in Table 2, above) that are provided sterile are packaged as follows:





The above listed materials are identical to those currently in use for Zimmer sterile instrument packages. Only the configuration of the interior components changes to correspond to the particular product being packaged.



The above listed materials are identical to those currently in use for other similar sterile instrument packages.

CONFIDENTIAL AND PROPRIETARY INFORMATION EXEMPT FROM DISCLOSURE UNDER FOI



FDA Question 4

	_	· · · · · · · · · · · · · · · · · · ·	Yes	N/A	No			
28.		he device, and/or accessory, and/or a component is provided						
I		'ile:	**					
		ect "N/A" if no part of the device, accessories, or						
	con	nponents is provided sterile, otherwise complete a-f below.						
	а.	Sterilization method is stated for each component	X					
		(including parameters such as dry time for steam	7					
		sterilization, radiation dose, etc.)		÷.				
	b.		X		Ì			
		A description of method to validate the sterilization cycle						
		(e.g., half-cycle method and full citation of FDA-						
		recognized standard, including date) is provided		_	x			
	C.	For devices sterilized using chemical sterilants such as			^			
i		ethylene oxide (EO) and hydrogen peroxide, submission						
	·	states maximum levels of sterilant residuals remaining on						
		the device and sterilant residual limits.						
		Select "N/A" if not sterilized using chemical sterilants.						
	d		X					
		Submission includes description of packaging and						
		packaging contents (e.g., if multiple devices are included						
		within the same package)						
	e.	Sterility Assurance Level (SAL) stated	X					
	f.				х			
		If device is blood-contacting, a permanent implant, or						
		contacts cerebrospinal fluid, or device is labeled "non-						
		pyrogenic," submission contains a description of the						
		endotoxin method used to make a determination (e.g.,						
		LAL), endotoxin release specification (e.g., 20 EU/device),						
		and a rationale for the specification.						
	1	Select "N/A" if device is not blood-contacting, not a						
		permanent implant, does not contact cerebrospinal fluid,						
		and is not labeled "non-pyrogenic." Select "N/A" if a						
		rationale for omission is provided.						
	Comments: (c) The sponsor has provided a promissory note for EO residuals.							
	Please provide EO residuals remaining on the device following sterilization							
		validation.						
	(f) The sponsor does not provide a rationale for omission of Pyrogen Testing.							
					-			



Zimmer Response to FDA Question 4

- (c) As discussed during the 15 January 2013 phone conversation between Joanna Surma (Zimmer, Inc.) and Peter Allen (FDA), the EO residual data will be provided to Peter Allen via email following each of the following testing phases:
 - EO Residuals after 1 Sterilization Cycle: expected by mid-February
 - EO Residuals after 2 Sterilization Cycles: expected by early March
 - EO Residuals after 3 Sterilization Cycles: expected by early- to mid-March
- (f) Page 87 of the premarket notification submission K123862 stated the following rationale for omission of Pyrogen Testing (additionally, Zimmer makes no "pyrogen-free" or "non-pyrogenic" label claims):

"This device is not labeled as nonpyrogenic. Per USP XXIII (161), requirements for specified endotoxin levels do not apply to orthopaedic implants."



FDA Question 5

		Yes	N/A	No
32.	Submission identifies contact classification (e.g., surface- contacting, less than 24 hour duration)			x
	Comments: Sponsor does not state the contact classificati	on for	the	
	instruments.			

Zimmer Response to FDA Question 5

Of the Class II Accessories to the Implant listed in Section 11a. (Device Description) of the premarket notification submission K123862, all are *external communicating devices* contacting *tissue/bone* with *limited exposure* (less than 24 hour duration), according to definitions presented in Section 5 of ISO 10993-1.



P.O. Box 708 Warsaw, IN 46581-0708 574 267-6131

510(k) Summary of Safety and Effectiveness

S	ponsor:
-	ponour.

Zimmer, Inc. P.O. Box 708 Warsaw, IN 46581-0708

Joanna L. Surma Associate Project Manager, Regulatory Affairs Telephone: (574) 371-1642 Fax: (574) 372-4605

Date:

Trade Name:

Contact Person:

Common Name:

Classification Name and Reference:

Product Code:

Predicate Device:

Device Description:

13 December 2012

Zimmer Nexel Total Elbow

Total Elbow Prosthesis

Prosthesis, Elbow, Constrained, Cemented 21 CFR § 888.3150

JDC

Coonrad/Morrey Total Elbow, manufactured by Zimmer, K001989, cleared 25 July 2000

Coonrad/Morrey Total Elbow, manufactured by Zimmer, K053189, cleared 9 December 2005

The Zimmer Nexel total Elbow is a total elbow prosthesis designed for use with bone cement. It is available in multiple sizes and in right and left configurations.

How the Device Functions: The Zimmer Nexel Total Elbow is an implant designed to replace the articulating surfaces of and restore motion to the humeroulnar (elbow) joint. The implant is a semi-constrained device assembly and consists of the following components: humeral component, ulnar component, humeral bearing-A, ulnar bearings-B, axle pin and humeral screws.

The Following Scientific Concepts, Design Features and Physical Properties form the Basis for the Zimmer Nexel Total Elbow: The humeral component has a humeral stem designed to be implanted with bone cement into the patient's humeral medullary canal, an anterior flange designed to accept a

510(k) Summary of Safety and Effectiveness - Zimmer Nexel Total Elbow

bone graft and limit torsional and posterior migration, a humeral yoke with rounded corners to avoid the creation of stress risers within the medial and lateral humeral supracondylar columns, and plasma spray region to enhance fixation to bone cement within the medullary canal, and to improve fatigue strength. The ulnar component has an ulnar stem designed to be implanted with bone cement into the patient's ulnar medullary canal, an ulnar eye that is both highly polished and nitrogen-enriched to limit wear of the apposing polymer bearings, and plasma spray region to enhance fixation to bone cement within the medullary canal. Bearings A and B are designed to broadly distribute joint reaction forces.

Materials Used: The humeral and ulnar components are made of a titanium alloy, the bearings A and B are made of Vitamin E highly cross-linked polyethylene, and the humeral screws and axle pin are made of a cobalt-chromium-molybdenum alloy.

Comparison to the Predicate: The proposed device (Zimmer Nexel Total Elbow) and the predicate device (Coonrad/Morrey Total Elbow) have the same intended use and similar indications for use. The proposed device humeral and ulnar components are very similar in terms of materials used and design/dimensions to the predicate device humeral and ulnar components. The primary differences between the proposed and predicate devices are that the proposed device bearing components are made of Vitamin E highly cross-linked polyethylene (the predicate device bearing components are made of ultra-high molecular-weight polyethylene), and the proposed device has bearings that articulate on both the outer and inner diameters of the ulnar eye . (the predicate device has bearings that articulate on the inner diameter of the ulnar eye). Additionally, the proposed device uses screws and a 1-part axle pin to fix the bearings in place, while the predicate device uses a 2-part (snap-fit) axle pin.

Intended Use:

Indications for use include:

- Elbow joint destruction which significantly compromises the activities of daily living
- Post-traumatic lesions or bone loss contributing to elbow instability
- Ankylosed joints, especially in cases of bilateral ankylosis from causes other than active sepsis
- Advanced rheumatoid, post-traumatic, or degenerative arthritis with incapacitating pain
- Instability or loss of motion when the degree of joint or soft tissue damage precludes reliable osteosynthesis
- Acute comminuted articular fracture of the elbow joint surfaces that precludes less radical procedures, including 13-C3 fractures of the distal humerus
- Revision arthroplasty

Caution: This device is intended for cemented use only.

510(k) Summary of Safety and Effectiveness - Zimmer Nexel Total Elbow

Comparison to Predicate Device:	The Zimmer Nexel Total Elbow is substantially equivalent to the predicate device in terms of form and function. The Zimmer Nexel Total Elbow and the predicate device share similar intended uses and indications for use.
Performance Data:	 Non-Clinical Performance Testing Conducted: Stem Fatigue Testing Wear Testing Durability Testing Modular Connection Fatigue Testing
	Non-Clinical Performance Testing Conclusions:
	Non-clinical testing demonstrated that the New Zimmer Total

Non-clinical testing demonstrated that the New Zimmer Total Elbow meets performance requirements as defined by Design Control activities and is substantially equivalent to the predicate device in terms of safety and efficacy.

In this case, clinical data and conclusions were not needed to demonstrate substantial equivalence.

Records processed under FOIA request #2016-4653; Released by 20RH on 07/16/2017



Traditional 510(k) Premarket Notification

Predicate Device(s)

The predicate and design basis for the Zimmer Nexel Total Elbow (the proposed device) is the Coonrad/Morrey Total Elbow, manufactured by Zimmer. Clearance letters (K001989 and K053189 cleared 25 July 2000 and 9 December 2005, respectively) for the predicate device are included at the end of this section.

Similarities and differences between the proposed implants and the predicate implants are tabulated below.

Substantial Equivalence Comparison

Indications for Use

	Proposed Device: The Zimmer Nexel Foud Elbow	Predicate Device: The Coonrad/Morrey Total Elbow (K001918, K053189)
Intended Use	To be used in conjunction with bone cement to replace the articulating surfaces of the humeroulnar joint.	To be used in conjunction with bone cement to replace the articulating surfaces of the humeroulnar joint.
Indications for Use	 Elbow joint destruction which significantly compromises the activities of daily living Post-traumatic lesions or bone loss contributing to elbow instability Ankylosed joints, especially in cases of bilateral ankylosis from causes other than active sepsis Advanced rheumatoid, post-traumatic, or degenerative arthritis with incapacitating pain Instability or loss of motion when the degree of joint or soft tissue damage precludes reliable osteosynthesis Acute comminuted articular fracture of the elbow joint surfaces that precludes less radical procedures, including 13-C3 fractures of the distal humerus Revision arthroplasty CAUTION: This device is intended for cemented use only. 	 Elbow joint destruction which significantly compromises the activities of daily living Post-traumatic lesions or bone loss contributing to elbow instability Ankylosed joints, especially in cases of bilateral ankylosis from causes other than sepsis Advanced rheumatoid or degenerative arthritis with incapacitating pain Instability or loss of motion when the degree of joint damage precludes less radical procedures Revision arthroplasty CAUTION: This device is intended for cemented use only.
Product Code	JDC	JDC

CONFIDENTIAL AND PROPRIETARY INFORMATION EXEMPT FROM DISCLOSURE UNDER FOI




ì

Traditional 510(k) Premarket Notification

Technology

	Feature	Proposed Device: The Zimmer Nevel Total Elbow	Predicate Device: The Coonrad/Morrey Total Elbow (K001918, K053189)
	Material	Wrought Tivanium Alloy (Ti-6Al-4V) with Titanium Plasma Spray	Wrought Tivanium Alloy (Ti-6Al-4V) with a porous coating of Titanium Beads
	Fixation Method	Intended for cemented use only	Intended for cemented use only
- General	Comparable Sizes	Size 4, Standard Flange, 100mm $\rightarrow \rightarrow \rightarrow$ Size 4, Standard Flange, 150mm $\rightarrow \rightarrow \rightarrow \rightarrow$ Size 5, Standard Flange, 100mm $\rightarrow \rightarrow \rightarrow$ Size 5, Standard Flange, 150mm $\rightarrow \rightarrow \rightarrow$ Size 5, Long Flange, 150mm $\rightarrow \rightarrow \rightarrow \rightarrow$ Size 6, Standard Flange, 100mm $\rightarrow \rightarrow \rightarrow$ Size 6, Standard Flange, 100mm $\rightarrow \rightarrow \rightarrow$ Size 6, Standard Flange, 150mm $\rightarrow \rightarrow \rightarrow$ Size 6, Long Flange, 150mm $\rightarrow \rightarrow \rightarrow \rightarrow$ Size 6, Long Flange, 200mm $\rightarrow \rightarrow \rightarrow \rightarrow \rightarrow$	 → Size Extra Small, Standard Flange, 4in → Size Extra Small, Standard Flange, 6in → Size Small, Standard Flange, 4in → Size Small, Standard Flange, 6in → Size Small, Long Flange, 6in → Size Regular, Standard Flange, 4in → Size Regular, Standard Flange, 6in → Size Regular, Long Flange, 6in → Size Regular, Long Flange, 6in → Size Regular, Long Flange, 6in
nponent	Sterility	Provided Sterile via Gamma Irradiation	Provided Sterile via Gamma Irradiation
Humeral Component – General	Design Features	 Humeral stem designed to be implanted with bone cement into the patient's humeral medullary canal Anterior flange designed to accept a bone graft and limit torsional and posterior migration Humeral yoke with rounded corners to avoid the creation of stress risers within the medial and lateral humeral supracondylar columns Plasma spray region to enhance fixation to bone cement within the medullary canal, and to improve fatigue strength Features to accept Humeral Screws, Axle Pin and polymer Bearings A and B 	 Humeral stem designed to be implanted with bone cement into the patient's humeral medullary canal Anterior flange designed to accept a bone graft and limit torsional and posterior migration Humeral yoke with squared corners Titanium Bead coating region to enhance fixation to bone cement within the medullary canal Features to accept Axle Pin and polymer Bearings






CONFIDENTIAL AND PROPRIETARY INFORMATION EXEMPT FROM DISCLOSURE UNDER FOI



	Feature	Proposed Device: The Zimmer Nevel Total Fibow	Predicate Deviçe: The Coonrad/Morrey Total Elbow (K001918, K053189)
	Material	Wrought Tivanium Alloy (Ti-6Al-4V) with Titanium Plasma Spray	Wrought Tivanium Alloy (Ti-6Al-4V) with Titanium Plasma Spray
	Fixation Method	Intended for cemented use only	Intended for cemented use only
Ulnar Component – General	Left, Size 4, 115mm $\rightarrow \rightarrow $	 → Left, Size Extra Small, Standard → Left, Size Extra Small, Long → Left, Size Small, Standard → Left, Size Small, Long → Left, Size Regular, Standard → Left, Size Regular, Long → Right, Size Extra Small, Standard → Right, Size Extra Small, Long → Right, Size Small, Standard → Right, Size Small, Long → Right, Size Small, Long → Right, Size Regular, Standard 	
Design	Sterility	Provided Sterile via Gamma Irradiation	Provided Sterile via Gamma Irradiation
	Design Features	 Left and right side-specific Ulnar stem designed to be implanted with bone cement into the patient's ulnar medullary canal Ulnar eye that is both highly polished and nitrogen-enriched to limit wear of the apposing polymer bearings Plasma spray region to enhance fixation to bone cement within the medullary canal 	 Left and right side-specific Ulnar stem designed to be implanted with bone cement into the patient's ulnar medullary canal Ulnar eye that is highly polished to limit wear of the apposing polymer bearings Plasma spray region to enhance fixation to bone cement within the medullary canal

CONFIDENTIAL AND PROPRIETARY INFORMATION EXEMPT FROM DISCLOSURE UNDER FOI





CONFIDENTIAL AND PROPRIETARY INFORMATION EXEMPT FROM DISCLOSURE UNDER FOI

23




	Feature	Proposed Device: The Zimmer Nexel Fotal Etbow	Predicate Device: The Cooncad/Morrey Total Elbow (K001918, K053189)
	Material	Vitamin E HXPE (Vitamin E Highly Cross-Linked Polyethylene)	UHMWPE (Ultra-High Molecular-Weight Polyethylene)
. _	Fixation Method	Press-Fit to Humeral Component	Sliding-Fit to Humeral Component Press-Fit to Ulnar Component
nera	Sterility	Provided Sterile via Ethylene Oxide	Provided Sterile via Gamma Irradiation
Polyethylene Components - General	Design Features	 Bearings A and B: Designed to broadly distribute joint reaction forces Designed to articulate against the Eye of the Ulnar Component Features to allow for press-fit mate to Humeral Component Bearings B only: Features to allow for press-fit mate to Axle Pin Designed to articulate against the inner-diameter and the medial and lateral surfaces of the Eye of the Ulnar Component Bearing A only: Designed to articulate against the outer-diameter of the Eye of the Ulnar Component 	 Humeral Bushings: Designed to articulate with the medial and lateral sides of the Eye of the Ulnar Component Features to allow for sliding-fit mate to the Humeral Component and sliding-fit mate to the Axlc Pin (also known as Hinge Pin, also known as Snap Pin) Ulnar Bushing: Features to allow for a press-fit mate with the Eye of the Ulnar Component Designed to articulate against the Humeral Bushings

CONFIDENTIAL AND PROPRIETARY INFORMATION EXEMPT FROM DISCLOSURE UNDER FOI

ł







	Feature	Proposed Device: The Zimmer Nevel Total Efform	Predicate Device: The Coonrad/Morrey Total Elbow (K001948, K053189)	
	Material	Zimaloy (Cobalt-Chromium-Molybdenum Alloy)	Zimaloy (Cobalt-Chromium-Molybdenum Alloy)	
Axle Pin	Fixation Method	 Press-Fit to Bearings B Compression fit between the taper of the Humeral Screws and the yoke of the Humeral Component 	 Sliding-Fit to Bearings Snap-Fit between inner and outer Axle Pin components 	
	Sizes	Size 4 (b) ' Long ' Diameter • Solid Material Throughout Size 5/6 (b)(4) Long Diameter • Solid Material Throughout	Size Extra Small (b)(4) Long Outer Diameter Wall Thickness Size Small/Regular (b) (4) Outer Diameter Wall Thickness	
	Sterility	Provided Sterile via Ethylene Oxide	Provided Sterile via Gamma Irradiation	
	Design Features	 Designed to remain in a fixed position (not translate nor rotate) throughout the entire humeroulnar joint range of motion Machined complete as a single, solid component 	 Designed to freely rotate and slide Machined as 2 (an inner and an outer) hollow components that snap together intraoperatively 	
	Material	Zimaloy (Cobalt-Chromium-Molybdenum Alloy)		
	Fixation Method	Threads engage with the Spiralock® threads of the Humeral Component		
	Sizes	Single size mates with all humeral component sizes and both axle pin sizes	Not Applicable: No comparable part in the predicate system	
crews	Sterility	Provided Sterile via Gamma Irradiation		
Humeral Screws	Design Features	 Designed to remain in a fixed position throughout the entire humeroulnar joint range of motion (not directly loaded by any compressive joint reaction force during flexion or extension of the humeroulnar joint) Headless design intended to ensure the constant application of a clamp-load to the axle-pin and to avoid soft tissue impingement or irritation associated with a protruding screw head Non-threaded taper to (along with the Humeral Component) hold Axle Pin in compression Pilot feature to avoid cross-threading 		

CONFIDENTIAL AND PROPRIETARY INFORMATION EXEMPT FROM DISCLOSURE UNDER FOI

1



		iemi-Constraint Characteristics	Proposed Device: The Zimmer Nevel Total Elbow	Predicate Device: The Coonrad/Morrey Fotal Elbow (K081918, K053189)	
nponent)		Varus/Valgus Angular Laxity	Neutral to Maximum Varus: (b) Neutral to Maximum Valgus: (b) Total Included Angle: (b)	Neutral to Maximum Varus: (b)(Neutral to Maximum Valgus: (b)(Total Included Angle: (b)	
to Humeral Con		Internal/External Rotational Laxity	Neutral to Maximum Internal Rotation: (b) Neutral to Maximum External Rotation:	Neutral to Maximum Internal Rotation: (b)(Neutral to Maximum External Rotation: (b)(4 Total Included Angle: (b)(-
Motion respect	Pistoning	Proximal/Distal Translational Laxity	Neutral to Maximum Proximal (b)(4) Neutral to Maximum Distal: Total Proximal/Distal Laxity:	 Neutral to Maximum Proximal (b)(4) Neutral to Maximum Distal: Total Proximal/Distal Laxity: 	Greated allow chat of tron to a weat
ange of lent with	Pisto	Anterior/Posterior Translation Laxity	Neutral to Maximum Anterior: Neutral to Maximum Posterior Total Anterior/Posterior Laxity	Neutral to Maximum Anterior: Neutral to Maximum Posterior Total Anterior/Posterior Laxity	concern, and
rained R t movem		Medial/Lateral Translation Laxity	Neutral to Maximum Medial: Neutral to Maximum Lateral: Total Medial/Lateral Laxity:	Neutral to Maximum Medial: Neutral to Maximum Lateral: Total Medial/Lateral Laxity:	This could ssed was could ssed in loomet-
Semi-Constrained Range of Motion (defined in terms of Ulnar Component movement with respect to Humeral Component)		Flexion/Extension Range defined in terms of degrees flexion)	All Humeral and Ulnar combinations: (b)(4)	Humeral size XS with Ulnar size XS: (b)(4) Humeral size XS with Ulnar size S: (b)(4) Humeral size XS with Ulnar size REG: (b)(4) Humeral size S with Ulnar size XS: (b)(4) Humeral size S with Ulnar size REG: (b)(4) Humeral size REG with Ulnar size XS: (b)(4) Humeral size REG with Ulnar size S: (b)(4) Humeral size REG with Ulnar size S: (b)(4) Humeral size REG with Ulnar size REG: (b)(4) Humeral size REG with Ulnar size REG with Ulnar size REG: (b)(4) Humeral size REG wit	~**
			لر	Is this a concern. constraint, stability, loutening),
EXEM	PT	FROM DISCLOSUR	E UNDER FOI	constraint, stability, 10020 270 Stats, tracture] TUS@fda.bbs.gov.or.301.706.8118	



Why Any Differences Do Not Render the Device NSE

The major technological differences between the proposed device (Zimmer Nexel Total Elbow) and the predicate device (the Zimmer Coonrad/Morrey Total Elbow) are listed below along with the performance tests conducted to <u>demonstrate substantial equivalence</u> (the listed performance tests and their conclusions are described in more detail in *Section 18a Performance Testing – Bench* of this premarket notification submission):

Device Component	Major Difference Between Proposed and Predicate Devices	Performance Test Conducted to Demonstrate Substantial Equivalence
Humeral Component	 Distal design features for implant assembly Plasma Spray 	 Durability Testing Modular Connection Fatigue Testing Fatigue Testing
Ulnar Component	- Nitrogen enrichment	- Wear Testing - Fatigue Testing
Polyethylene Components	 Polyethylene material Surfaces against which components are intended to articulate 	 Wear Testing Durability Testing
Axle Pin	- Single Solid vs. 2-Part Hollow construction	 Wear Testing Durability Testing Modular Connection Fatigue Testing
Humeral Screws	- Existence of screws	 Wear Testing Durability Testing Modular Connection Fatigue Testing

Current questions of safety and effectiveness (clinically observed failure modes) associated with elbow prostheses are humeral and ulnar component loosening, humeral and ulnar component fracture, polyethylene wear, system (implant assembly) durability and system (implant assembly) disassociation (also known as modular connection fatigue). To demonstrate that the proposed device is at least as safe and effective as the predicate device despite the compounding effects of these technological differences and previously outlined dimensional differences, fatigue (to assess resistance to fracture), wear, durability and modular connection fatigue performance tests were conducted (testing conclusions are described in more detail in *Section 18a Performance*

CONFIDENTIAL AND PROPRIETARY INFORMATION EXEMPT FROM DISCLOSURE UNDER FOI



A .- 1

4:54

Are there other

4 (65°

Elexin?

predictors

Testing – Bench of this premarket notification submission). Stem loosening was not the focus of a performance test for the New Zimmer Total Elbow, because the variables (cementing technique, stem geometry, stem material and stem surface finish) associated with stem loosening are consistent from the predicate to the proposed device; and, therefore the proposed device (the New Zimmer Total Elbow) is expected to be as safe and effective as the predicate device (the Coord/Morrey Total Elbow) in terms of ability to resist stem loosening.

As reflected in the dFMEA documentation associated with the proposed device, the Zimmer Nexel Total Elbow does not raise any new/different types of questions of safety or effectiveness.

The proposed and predicate devices have the exact same intended use. The proposed and predicate devices share similar indications for use; where, the only changes in wording of the indications for use are that the proposed device has indications for use that are intended to add specificity to the indications for use of the proposed device. This added specificity does not change the intended meaning of the indications for use of the proposed device as compared to those of the predicate.

Because the aforementioned differences between the proposed and predicate devices do not constitute a new intended use, affect safety or effectiveness or raise different questions of safety and effectiveness, these differences *do not* render the proposed device NSE (not substantially equivalent).



Performance Specifications

	Test Results		
Test	Proposed Device: The Zimmer Nexet Fotal Fibow	Predicate Device: The Coonrad/Morrey Total Elbow (K001918, K053189)	
Fatigue Testing: Humeral Stem (Plasma Spray Region)	Five Zimmer Nexel Total Elbow humeral stems achieved (b)(4) load cycles at (b)(6) without fracture.	The predicate device humeral component was fatigue tested at various load levels, and the fatigue strength at (b)(d) cycles was determined to be (b)(
Fatigue Testing: Ulnar Stem (Mid-Stem Region)	Five Zimmer Nexel Total Elbow ulnar stems achieved (b)(4) load cycles in the mid-stem region at a load of (b) without fracture.	The predicate device ulnar component was fatigue tested at various load levels, and the fatigue strength at $(b)(4)$ load cycles in the mid-stem region was determined to be $7(5)7$	
Fatigue Testing: Ulnar Stem (Plasma Spray Region)	Five Zimmer Nexel Total Elbow ulnar stems achieved (b)(4) load cycles in the plasma spray region at a load of (b)(without fracture.	The predicate device ulnar component was fatigue tested at various load tevels, and the fatigue strength at $(b)(4)$ load cycles in the plasma spray region was determined to be (b)	
Wear Testing	The Zimmer Nexel Total Elbow mean gravimetric wear rate $(b)(4)$ at $(b)(4)$ load cycles) was 91% less than that of the predicate.	The predicate device mean gravimetric wear rate was (b) mg/Mc at (b)(4) load cycles.	
Durability Testing	The Zimmer Nexel Total Elbow device achieved (h)(4) run-outs at an equivalent of (b) f weight-in-hand.	The run out load for the predicate device was found to be equivalent to (6) (4) weight-in-hand.	
Modular Connection Fatigue Testing	The median fatigue strength of the non- articulating, mechanically locked, modular Zimmer Nexel Total Elbow implant components was no less than (b)(Not Applicable: The Zimmer Coonrad/Morrey Total Elbow does not have non-articulating, mechanically-locked, modular implant components comparable to those of the proposed device.	





Manual Orthopaedic Surgical Instruments: Recommendations for Care, Cleaning, Maintenance and Sterilization

Page intentionally left blank.

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017

Manual Orthopaedic Surgical Instruments

Table of Contents

1.	Purpose
2.	Scope
3.	Glossary
4.	Acronyms
5.	Symbols
6.	Considerations
7.	Processing Category Codes
8.	Instructions
Q.	
	A. Warning and Precautions
	C. Limitations and Restrictions on Reprocessing
	D. Point of Use Preparation for Reprocessing
	E. Preparation Before Cleaning
	F. Preparation of Cleaning Agents
	G. Manual Cleaning and Disinfection Instructions
	H. Combination Manual/Automated Cleaning and Disinfection Instructions
	I. Automated Cleaning and Disinfection Instructions
	J. Inspection, Maintenance, Testing and Lubrication
	K. Sterile Packaging
	L. Sterilization
	M. Storage Recommendations
9.	Hospital Responsibilities for Zimmer Loaned Instrument Sets
10	Customer Service Information
11.	References
	Table 1. Cleaning/Disinfection Options 9
	Table 2. Typical Automated Washer/Disinfector Cycle for Surgical Instruments 10
	Table 3. Recommended Steam Sterilization Parameters 12
	Appendix 1. Cleaning/Disinfection Processes
	Chart 1. Manual Cleaning/Disinfection Procedure15
	Chart 2. Combination Manual/Automated Cleaning/Disinfection Procedure15
	Chart 3. Typical Automated Washer/Disinfector Standard Instrument Cleaning/Disinfection Cycle 15



1. Purpose

These instructions are recommended for the care, cleaning, maintenance and sterillzation of reusable Zimmer orthopaedic manual surgical instruments. This document is intended to assist health care personnel in safe handling practices, effective reprocessing and maintenance of Zimmer reusable devices.

The manual is intended to assist the hospital and central supply management in developing procedures for safe and effective reprocessing of Zimmer instrument sets.

Hospital personnel, including those in receiving and central sterile supply departments (CSSD), as well as in the operating room (OR) may be directly involved in handling instruments purchased from Zimmer or on a loan basis as consignment instruments. Hospital directors and other management in each of these departments should be informed of these instructions and recommendations to ensure safe and effective reprocessing and to prevent damage or misuse of reusable devices.

2. Scope

This Instruction manual provides Information on the care, cleaning, disinfection, maintenance and sterilization of manual surgical instruments and is <u>applicable</u> to all reusable medical devices manufactured and/or distributed by Zimmer, Inc.

This information is also <u>applicable</u> to single-use medical devices manufactured by Zimmer that are supplied nonsterile but are intended to be used in a sterile state. These devices are single-use but can be reprocessed if <u>not used</u> (e.g. screws, plates, etc.). This also includes single-use devices packaged and sold sterile but removed from packaging and placed in kits.

Note: <u>not used</u> refers to those single-use components that have not been in contact with blood, bone, tissue or other body fluids. Any unused, single-use device that has been exposed to blood, bone, tissue or body fluids must not be reprocessed or resterilized and must be discarded. Devices that cannot be reused may be labeled with the following symbol:



Do not reuse

This information is <u>not applicable</u> to single-use devices that are sold sterile and cannot be resterilized (e.g. osteotome blades).

Devices that cannot be resterilized may be labeled with the following symbol:



Do not resterilize

This instruction manual is <u>not applicable</u> to air driven or electrically powered equipment. However, it is applicable to functional attachments (e.g. reamers and drill bits) that are connected to powered equipment for use. Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017

Manual Orthopaedic Surgical Instruments

3. Glossary

<u>Chemical</u>: a formulation of compounds intended for use in reprocessing.

Note: This includes detergents, surfactants, rinse aids, disinfectants, enzymatic cleaners and sterilants.

<u>Cleaning:</u> the removal of contamination from an item to the extent necessary for further processing.

<u>Contaminated</u>: State of having been actually or potentially in contact with microorganisms.

Decontamination: the use of physical or chemical means to remove, inactivate, or destroy bloodbome pathogens on a surface or item to the point where they are no longer capable of transmitting infectious particles and the surface or items is rendered safe for handling or disposal.

Disinfection: a process used to reduce the number of viable microorganisms on a product to a level previously specified as appropriate for its further handling or use.

Note: Cleaning and disinfection are often conducted in the same step.

Manual cleaning: cleaning without the use of an automated washer or washer/disinfector.

<u>Processing/Reprocessing</u>: activity including cleaning, disinfection and sterilization, necessary to prepare a new or used medical device for its intended use.

Sterile: free from all viable microorganisms.

<u>Sterilization</u>: a validated process used to render a device free from all forms of viable microorganisms.

Note: In a sterilization process, the nature of microbiological death is described by an exponential function. Therefore, the presence of microorganisms on any individual item may be expressed in terms of probability. While this probability may be reduced to a very low number, it can never be reduced to zero. This probability can only be assured for validated processes.

<u>Washer/Disinfector</u>: a machine intended to clean and disinfect medical devices and other articles used in the context of medical, dental, pharmaceutical, and veterinary practice.

4. Acronyms

- BI = biological indicator
- CJD = Creutzfeldt-Jakob Disease
- CSSD = Central Sterile Supply Department
- OR = operating room
- PPE = personal protective equipment
- SAL = sterility assurance level
- TSE = Transmissible Spongiform Encephalopathy

5. Symbols



180 16223 3.3

Do not reuse







Do not resterilize



Caution or Instructions for Use



6. Considerations

This instruction manual pertains to all Zimmer reusable surgical instruments and should be studied carefully. This manual supercedes Zimmer and Centerpulse Instrument manuals published prior to January 2006.

The user/processor should comply with local laws and ordinances in countries where reprocessing requirements are more stringent than those detailed in this manual.

New and used instruments must be thoroughly processed according to these instructions prior to use.

During musculoskeletal surgery, instruments become contaminated from blood, tissue, bone chips and marrow. The instruments may also become contaminated with body fluids containing hepatitis virus, HIV or other etiological agents and pathogens. All health care workers should become familiar with the necessary Universal Precautions of preventing injuries caused by sharp instruments when handling these devices during and after surgical procedures and during reprocessing.

It should be noted that saline and other irrigation fluids are often used in copious amounts during surgical procedures and will exert a corroding effect on instruments.

Orthopaedic surgery requires instruments which are heavy and have multiple components, articulating or rotating parts, removable handles, plastic replacement parts, and series of gauges or other measuring devices in graduated sizes. Devices are usually supplied in sets and subdivided into trays and cases in which the devices may be arranged by size or in the order needed for a Specific surgical procedure.

Hospitals must assume responsibility for cleaning, disinfection, packaging and sterilization of all loaner instrument sets before returning them to Zimmer. However, the next user must also inspect the set upon receipt to verify that instruments have, in fact, been adequately cleaned and decontaminated before repeating reprocessing procedures to prepare the loaner set for subsequent reuse. Zimmer cannot guarantee that sterility was attained by the previous user and has been maintained during transit. Zimmer representatives often open and inspect instrument sets between users, which will, of course, compromise sterility and require complete reprocessing prior to subsequent use. This manual includes instructions for all Zimmer reusable devices including legacy Centerpulse instruments marked with reprocessing category codes [a, a+, b, b+, c]. See Section 7 of this manual for further explanation of reprocessing codes. All Zimmer devices may be safely and efficiently reprocessed using the manual or combination manual/ automated cleaning instructions outlined in this manual.

Core orthopaedic instrument sets must be complete and in good condition to be used correctly. Optional devices may be available on request from your Zimmer representative. To maintain instruments properly it is important to consider the following information and processing instructions:

- Warnings and precautions
- Instrument set completeness and functionality
- Reprocessing limitations and/or restrictions
- Preparation for reprocessing at the point of use
- Preparation for cleaning (including assembly/disassembly as necessary)
- · Cleaning, disinfection and drying
- Maintenance, inspection, testing and lubrication
- Sterile packaging
- Sterilization
- Storage

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017

Manual Orthopaedic Surgical Instruments



7. Processing Category Codes

The following codes are etched on some instruments and provide information useful in the selection of cleaning agents with appropriate pH. Zimmer recommends that all reusable devices (regardless of etching) be processed in accordance with the manual or combination manual/automated cleaning instructions contained in this instruction manual.

a	Steel/metal instruments without cannulated bores/lumens or non-metal/polymer handles, or other components (e.g. retractors, drills, testing trays, rasps, scissors, clamps, exploring hooks, compression forceps, skin bridge elevators, guide wires, etc.). These devices are tolerant of alkaline cleaning agents when followed by acidic neutralization and thorough rinsing. These devices can be cleaned with rust-removal agents approved for surgical instruments in the presence of rust or corrosion.
a+	Steel/metal instruments with cannulated bores/lumens but without non-metal/polymer handles or other components (e.g. drills with elongated holes, belt tensioning pulleys, bone joint reamers, extractor cases). These devices are tolerant of alkaline cleaning agents when followed by acidic neutralization and thorough rinsing. These devices can be cleaned with rust-removal agents approved for surgical instruments in the presence of rust or corrosion. Cannulations and hollow spaces must be cleaned manually .
b	Instruments made of polymers or metal instruments paired with polymer components (e.g. testing trays for flat profiles, chisels with non-metal handles, awis, dissectors, femur dilatators, pyramidal chisels/rasps). These devices are tolerant of alkaline cleaning agents when followed by acidic neutralization and thorough rinsing.
b+	Instruments with cannulated bores, made of polymers or metal instruments paired with polymer components (e.g. tibial mallets, flex screwdrivers, tibial dilatators, etc.). These devices are tolerant of alkaline cleaning agents when followed by acidic neutralization and thorough rinsing. Cannulations and hollow spaces must be cleaned manually .
C	Instruments made of titanium or aluminum alloys and/or having assembly/disassembly or other reprocessing aids (e.g. torque spanners, tibial aiming devices, pad cutters, instrument cases, trays and sterilization containers). These devices should be cleaned using the manual or combination manual/ automated cleaning procedures provided in this manual. These devices should not be exposed to alkaline cleaning agents.



8. Processing Instructions

These processing instructions are intended to assist the hospital and central supply management in developing procedures to attain the above goals, both for hospitalowned and for loaned instrument sets. This information is based on Zimmer testing and experience, material science, as well as widely accepted recommendations of the following organizations:

- American National Standards Institute (ANSI)
- American Society for Testing and Materials (ASTM)
- Association for the Advancement of Medical Instrumentation (AAMI)
- Association of Operating Room Nurses (AORN)
- Centers for Disease Control (CDC)
- German Instrument Working Group (AKI) Arbeitskreis Instrumenten-Aufbereitung
- International Standards Organization (ISO)
- International Association of Healthcare Central Service Material Management (IAHCSMM)
- National Health Service (NHS)
- Robert Koch Institute (RKI)
- Swissmedic
- World Health Organization (WHO)

Note: These instructions describe the necessary processing steps that new and used instruments must undergo to attain sterlity.

A. Warnings and Precautions

- Universal Precautions should be observed by all hospital personnel that work with contaminated or potentially contaminated medical devices. Caution should be exercised when handling devices with sharp points or cutting edges.
- Personal Protective Equipment (PPE) should be worn when handling or working with contaminated or potentially contaminated materials, devices and equipment. PPE includes gown, mask, goggles or face shield, gloves and shoe covers.

- Metal brushes or scouring pads must not be used during manual cleaning procedures. These materials will damage the surface and finish of instruments. Soft-bristled, nylon brushes and pipe cleaners should be used.
- Cleaning agents with low foaming surfactants should be used during manual cleaning procedures to ensure that instruments are visible in the cleaning solution. Manual scrubbing with brushes should always be performed with the instrument below the surface of the cleaning solution to prevent formation of aerosols and splashing which may spread contaminants. Cleaning agents must be easily and completely rinsed from device surfaces to prevent accumulation of detergent residue.
- . Do not place heavy instruments on top of delicate devices.
- Do not allow contaminated devices to dry prior to reprocessing. All subsequent cleaning and sterilization steps are facilitated by not allowing blood, body fluid, bone and tissue debris, saline, or disinfectants to dry on used instruments.
- Saline and cleaning/disinfection agents containing aldehyde, mercury, active chlorine, chloride, bromine, bromide, lodine or iodide are corrosive and should not be used. Instruments must not be placed or soaked in Ringers Solution.
- Mineral oil or silicone lubricants should not be used because they: 1) coat microorganisms; 2) prevent direct contact of the surface with steam; and 3) are difficult to remove.
- Only devices manufactured and/or distributed by Zimmer should be included in Zimmer instrument trays and cases. These validated reprocessing instructions are not applicable to Zimmer trays and cases that include devices that are not manufactured and/or distributed by Zimmer,
- Descaling agents that include morpholine should not be used in steam sterilizers. These agents leave residue which can damage polymer instruments over time.



B. Receiving Inspection — Instrument set content and functionality verification

- Upon receipt in the hospital, instrument sets should be inspected for completeness. Inspect for thumb, wing, set, or other types of screws; screw-in or other detachable handles; and auxiliary exchangeable parts such as blades, right/left attachments or heads. Many organizing cases have shadow graphs, outlines, catalog numbers, and instrument names or sizes silk-screened or otherwise marked on the case or tray.
- Orthopaedic surgical procedures follow a precise order in which the instruments are used. Also, many instruments have dimensional features which govern bone resections, determine implant sizes, and measure intramedullary canal sizes, depth of drill holes, angles of tube/plate, acetabular cup placements, etc. Therefore, it is very important that all requested sizes of a specific instrument series are available (specific instruments are routinely omitted from instrument sets due to infrequent use unless requested by the user). Contact your Zimmer representative if requested instruments have been omitted but are required for surgery.
- Markings on instruments used for measuring anatomical dimensions must be legible. These may include gauge markings, angles, inner or outer diameters, length or depth calibrations, and right/left indications. Notify your Zimmer representative if scales and other markings are not legible.

C. Limitations and Restrictions

 Neutral pH enzymatic and cleaning agents are recommended and preferred for cleaning Zimmer reusable devices. Alkaline agents with pH of 12 or less may be used to cleanstainless steel and polymer instruments in countries where required by law or local ordinance; or where prion diseases such as Transmissible Spongiform Encephalopathy (TSE) and Creutzfeldt-Jakob Disease (CJD) are a concern. It is critical that alkaline cleaning agents are completely and thoroughly neutralized and rinsed from devices.

Note: Drill bits, reamers, rasps and other cutting devices should be carefully inspected after processing with alkeline detergents to ensure that cutting edges are fit for use. Note: It is important to select enzymatic solutions intended for breakdown of blood, body fluids and tissues. Some enzymatic solutions are specifically for breakdown of fecal matter or other organic contaminants and may not be suitable for use with orthopaedic instruments.

- Repeated processing, according to the instructions in this manual has minimal affect on Zimmer reusable manual instruments unless otherwise noted. End of life for stainless steel or other metal surgical instruments is normally determined by wear and damage due to the intended surgical use and not to reprocessing.
- Automated cleaning using a washer/disinfector alone may not be effective for orthopaedic instruments with lumens, cannulations, blind holes, mated surfaces and other complex features. A thorough, manual or combination manual/automated cleaning process is recommended.
- Where applicable, multi-component instruments should be disassembled for cleaning. Disassembly, where necessary, is generally self evident. Care must be taken to avoid losing small parts. If a part is lost, notify your Zimmer representative when the instrument set is returned.
- Instruments must be removed from metal or polymer trays for manual and/or automated cleaning procedures. Do not clean instruments while in polymer or metal trays. Instrument trays, cases and lids must be cleaned separately from instruments. Non-sterile, single-use plate and screw implants are an exception to this rule. Plates and screws may remain in the tray or caddy for reprocessing.
- Polymers used in Zimmer instrument sets can be sterilized using steam/moist heat. Polymer materials have a limited useful life. If polymer surfaces turn "chalky," show excessive surface damage (e.g. crazing or delamination), or if polymer devices show excessive distortion or are visibly warped, they should be replaced. Notify your Zimmer representative if polymer devices need to be replaced.
- Most currently available polymers will not withstand conditions in washer/sterilizers that operate at temperatures equal to or greater than 141°C/285°F, and use live-steam jets as cleaning features. Severe surface damage to polymer devices will occur under these conditions.



- Soaking in disinfectants may be a necessary step to control certain viruses. However, these agents may discolor or corrode instruments (household bleach contains or forms chlorine and chloride in solution and has a corrosive effect similar to sallne). Disinfectants containing glutaraldehyde, or other aldehydes, may denature protein based contaminants, causing them to harden and making them difficult to remove. Where possible, soaking in disinfectants should be avoided.
- Steam/moist heat is the recommended sterilization method for Zimmer instruments.
- Ethylene Oxide (EO), Gas Plasma Sterilization and dry heat sterilization methods are not recommended for sterilization of Zimmer reusable instruments.
- Instrument with removable polymer sleeves must be disassembled for sterilization (e.g. acetabular reamer shaft with sleeve, side cutters, etc.)
- During initial steam sterilization runs some formaldehyde from polyformaldehyde surfaces may vaporize and become noticeable. This should not cause concern. After a few sterilization cycles, the odor should be no longer evident.
- While ethylene oxide sterilization may prolong the service life of certain polymers (e.g. polysulfone), this method of sterilization is not recommended unless aeration times are provided in specific package inserts. Large polyformaldehyde items (DELRIN, CELCON) have been found to require excessive outgassing times (a minimum of five days at elevated temperatures in a mechanical aerator); therefore, gas sterilization for polyformaldehyde products is contraindicated.
- Titanium and titanium alloy devices are especially susceptible to discoloration from steam impurities and detergent residues which form multi-colored surface layers of oxide deposits. Upon repeated sterilization these oxide layers, while not harmful to the patient, may become so dark that they can obscure graduation marks, catalog and lot numbers, and other stamped or etched information. Acidic, anti-corrosion agents may be used to remove this discoloration. Avoid frequent use of these agents.
- Use of hard water should be avoided. Softened tap water may be used for InItial rinsing. Purified water should be used for final rinsing to eliminate mineral deposits on instruments. One or more of the following processes may be used to purify water: ultra-filter (UF), reverse-osmosis (RO), delonized (DI), or equivalent.

D. Point of Use Preparation for Reprocessing

 Remove excess body fluids and tissue from instruments with a disposable, non-shedding wipe. Place instruments in a basin of distilled water or in a tray covered with damp towels. Do not allow saline, blood, body fluids, tissue, bone fragments or other organic debris to dry on instruments prior to cleaning.

Note: Soaking in proteolytic enzyme solutions facilitates cleaning, especially in instruments with complex features and hard-to-reach areas (e.g. cannulated and tubular designs, etc.). These enzymatic solutions break down protein matter and prevent blood and protein based materials from drying on instruments. Manufacturer's instructions for preparation and use of these solutions should be explicitly followed.

- Instruments should be cleaned within 30 minutes of use to minimize the potential for drying prior to cleaning.
- Used instruments must be transported to the central supply in closed or covered containers to prevent unnecessary contamination risk.

E. Preparation Before Cleaning

- Symbols or specific instructions etched on instruments or instrument trays and cases should be strictly followed.
- Where applicable, multi-component instruments should be disassembled for appropriate cleaning. Care should be exercised to avoid losing small screws and components. If a part is lost, notify your Zimmer representative when the instrument set is returned.
- Published instructions for use and surgical techniques and/or procedures may provide a supplemental source to illustrate assembly/disassembly instructions for specific Zimmer instruments.



F. Preparation of Cleaning Agents

- Neutral pH enzymatic and cleaning agents with low foaming surfactants are preferred and recommended by Zimmer.
 Alkaline agents with pH of 12 or less may be used in countries where required by law or local ordinance.
 Alkaline agents should be followed with a neutralizer and thorough rinsing.
- All cleaning agents should be prepared at the use-dilution and temperature recommended by the manufacturer.
 Softened tap water may be used to prepare cleaning agents. Use of recommended temperatures is important for optimal performance of cleaning agents.
- Dry powdered cleaning agents should be completely dissolved prior to use to avoid staining or corrosion of instruments.
- Fresh cleaning solutions should be prepared when existing solutions become grossly contaminated (bloody and/or turbid).

Method	Description	Section
Manual	Enzymatic soak and scrub followed by sonication	G
Combination Manual/Automated	Enzymatic soak and scrub followed by automated washer/ disinfector cycle	H
Automated (Washer/Disinfector)	Washer/disinfector cycle - not recommended for use without manual precleaning	1

Table 1. Cleaning/Disinfection Options

G. Manual Cleaning/Disinfection Instructions

- Completely submerge instruments in enzyme solution and allow to soak for 20 minutes. Use a soft-bristled, nylon brush to gently scrub the device until all visible soll has been removed. Particular attention must be given to crevices, lumens, mated surfaces, connectors and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft-bristled brush (i.e. pipe cleaner brush).
- Remove the device from the enzyme solution and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, holes and other difficult-toreach areas.
- Place prepared cleaning agents in a sonication unit. Completely submerge device in cleaning solution and sonicate for 10 minutes at 45-50kHz.
- 4. Rinse instrument in purified water for at least 3 minutes or until there is no sign of blood or soil on the device or in the rinse stream. Thoroughly and aggressively flush lumens, holes and other difficult-to-reach areas.
- 5. Repeat the sonication and rinse steps above.
- 6. Remove excess moisture from the instrument with a clean, absorbent and non-shedding wipe.

Note: If stainless steal instruments are stained or corroded, an acidic, anti-corrosion agent in an ultrasonic cleaner may be sufficient to remove surface deposits. Care must be taken to thoroughly rinse acid from devices. Acidic, anti-corrosion agents should only be used on an as needed basis.



H. Combination Manual/Automated Cleaning and Disinfection Instructions

 Completely submerge the instruments in enzyme solution and allow to soak for 10 minutes. Use a soft nylon-bristled brush to gently scrub the device until all visible soll has been removed. Particular attention must be given to crevices, lumens, mated surfaces, connectors and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft nylon-bristled brush (i.e. pipe cleaner).

Note: Use of a sonicator at 45-50kHz will aid in thorough cleaning of devices.

Note: Use of a syringe or water jet will improve flushing of difficult to reach areas and closely mated surfaces.

- Remove devices from the enzyme solution and rinse in purified water for a minimum of 1 minute. Thoroughly and aggressively flush lumens, blind holes and other difficult-toreach areas.
- Place instruments in a suitable washer/disinfector basket and process through a standard instrument washer/ disinfector cleaning cycle. The following minimum parameters are essential for thorough cleaning and disinfection.

Table 2. Typical Automated Washer/Disinfector Cycle for Surgical Instruments

Step	Description
1	2 minute prewash with cold tap water
2	20 second enzyme spray with hot tap water
3	1 minute enzyme soak
4	15 second cold tap water rinse (X2)
5	2 minutes detergent wash with hot tap water (64·66°C/146-150°F)
6	15 second hot tap water rinse
7	2 minute thermal rinse (80-93°C/176-200°F)
8	10 second purified water drise with optional lubricant (64-66°C/146-150°F)
9	7 to 30 minute hot air dry (116°C/240°F)

I. Automated Cleaning/ Disinfection Instructions

- Automated washer/disinfector systems are not recommended as the sole cleaning method for surgical instruments. Orthopaedic instruments should be cleaned following the manual or combination manual/automated cleaning procedure outlined in this manual except where specifically indicated,
- An automated washer/disinfector may be used as follow-up to the manual cleaning procedure above but is not required.
- 3. Simple instruments without multiple components, lumens/ cannulations, blind holes, mated surfaces, connectors and internal mechanisms or other complex features may be successfully cleaned and disinfected using a typical washer/disinfector cycle for surgical instruments as outlined in Table 2 of this manual. Devices should be thoroughly inspected prior to sterilization to ensure effective cleaning.

J. Inspection, Maintenance, Testing and Lubrication

- Carefully inspect each device to ensure that all visible contamination has been removed. If contamination is noted repeat the cleaning/disinfection process.
- 2. Visually inspect for completeness, damage and/or excessive wear.

Note: If damage or wear is noted that may compromise the function of the instrument, contact your Zimmer representative for a replacement.

- Check the action of moving parts (e.g. hinges, box-locks, connectors, sliding parts, etc.) to ensure smooth operation throughout the Intended range of motion.
- 4. Hinged, rotating, or articulating Instruments should be lubricated with a water soluble product (e.g. Instrument Milk or equivalent lubricant) intended for surgical instruments that must be sterilized. Some water-based instrument lubricants contain bacteriostatic agents which are beneficial. To remain effective, the expiration date specified by the manufacturer should be adhered to for both stock and use-dilution concentrations.

Note: Mineral oll or silicone lubricants should not be used because they 1) coat microorganisms; 2) prevent direct contact of the surface with steam; and 3) are difficult to remove.



Note: These lubrication instructions are not applicable to air-powered or electrical instruments. These devices have different requirements and should be lubricated according to the manufacturer's instructions.

- Check instruments with long slender features (particularly rotating instruments) for distortion.
- Where instruments form part of a larger assembly, check that the devices assemble readily with mating components.

K. Sterile Packaging

Packaging individual instruments

- Commercially available, medical grade steam sterilization pouches (e.g. paper, Tyvek[™] or equivalent) of the appropriate sizes may be used to double package single instruments. Ensure that the inner pouch Is large enough to contain the instrument without stressing the seals or tearing the packaging but small enough to be placed in a secondary pouch without compromising the integrity of the total package.
- Standard medical grade, steam sterilization wrap may be used to package individual instruments. The package should be prepared using the AAMI double wrap or equivalent method.

Note: If sterilization wraps are used they must be free of detergent residues. Reusable wraps are not recommended.

Packaging instrument sets in rigid trays and cases with lids

Safety Precaution: The total weight of a wrapped instrument tray or case should not exceed 11.4kg/25lbs. When placed in a sterilization container with gasketed lid the total package should not exceed 16kg/35lbs.

- Trays and cases with lids may be wrapped in standard medical grade, steam sterilization wrap using the AAMI double wrap method or equivalent.
- Trays and cases with lids may also be placed in an approved sterilization container with a gasketed lid for sterilization.

Note: Follow the sterilization container manufacturer's instructions for inserting and replacing sterilization filters in sterilization containers,

Instrument trays and cases with defined, preconfigured layouts

- Areas designated for specific devices shall contain only devices specifically intended for these areas.
- Optional ZImmer instruments should not be added to a preconfigured instrument tray or case unless a dedicated universal space or compartment has been included in the design and the guidelines described below for trays and cases without defined layouts or universal spaces can be applied.
- Only devices manufactured and/or distributed by Zimmer should be included in Zimmer instrument trays. These validated reprocessing instructions are not applicable to Zimmer trays that include devices that are not manufactured and/or distributed by Zimmer.

Universal instrument trays and cases without defined, preconfigured layouts or containing undefined universal spaces or compartments should only be used under the following conditions

- The total weight of a wrapped instrument tray or case should not exceed 11.4kg/25lbs. When placed in a sterilization container with gasketed lid the total sterilization package should not exceed 16kg/35lbs.
- Any device capable of disassembly must be disassembled prior to placement in the case.
- All devices must be arranged to ensure steam penetration to all instrument surfaces. Instruments should not be stacked or placed in close contact.
- The user must ensure that the instrument case is not tipped or the contents shifted once the devices are arranged in the case. Silicon mats may be used to keep devices in place.
- Only devices manufactured and/or distributed by Zimmer should be included in Zimmer instrument trays. Zimmer validated reprocessing instructions are not applicable to Zimmer trays that include devices that are not manufactured and/or distributed by Zimmer.



L. Sterilization Instructions

- See Table 3 for recommended minimum sterilization parameters that have been validated by Zimmer to provide a 10⁴ sterility assurance level (SAL).
- The hospital is responsible for in-house procedures for the reassembly, inspection, and packaging of the instruments after they are thoroughly cleaned in a manner that will ensure steam sterilant penetration and adequate drying. Provisions for protection of any sharp or potentially dangerous areas of the instruments should also be recommended by the hospital.
- Moist heat/steam sterilization is the preferred and recommended method for Zimmer orthopaedic instrument sets.

- Sterilizer manufacturer recommendations should always be followed. When sterilizing multiple instrument sets in one sterilization cycle, ensure that the manufacturer's maximum load is not exceeded.
- Instrument sets should be properly prepared and packaged In trays and/or cases that will allow steam to penetrate and make direct contract with all surfaces.
- Ethylene oxide or gas plasma sterilization methods should not be used unless package inserts for the applicable product specifically provide instructions for sterilization using these methods.
- Gravity displacement sterilization cycles are not recommended because cycle times are too long to be practical.

Cycle Type	Minimum	*Pressure	* Min Exposi	¹¹ Minimum	
-,,	Temperature		^{L1} Wrapped	"Unwnapped	Dry Time
^{1,1} UK Prevacuum/ Pulsating Vacuum	134°C 273°F	3bar 28.5psi	3 min	3 min	
^{1,1} Prevacuum/ Pulsating Vacuum	132ºC 270 °F	1.86bar 27psi	4 min	4 min	
** Prevacuum/ Pulsating Vacuum	134°C 273°F	3bar 28.5psi	18 min	18 min	30 minutes
<pre>Prevacuum/ Pulsating Vacuum</pre>	132°C 270°F	1.86bar 27psi	8 min	8 min	
¹² Gravity/Gravity Displacement	Not recommer are not practic		excessively lon	g sterilization o	ycles which

Table 3. Recommended Steam Sterilization Parameters

1 Minimum validated steam sterilization time required to achieve a 10*

- sterility assurance level (SAL).
 2 Minimum validated scenn sterifization <u>temperature</u> required to achieve a
- 10⁴ steniity assurance level (SAL). 3 Local or national specifications should be followed where steam steniization requirements are stricter or more conservative than those listed in this table.
- 4 Disinfection/steam sterilization parameters recommended by the World Health Organization (WHO) for reprocessing instruments where there is recommended by the standard st
- concern regarding TSE/CJÜ contamination. 5 For Universal Instrument Cases without defined load configurations.
- 6 Sea level
- 7 AAMI/AORN steam sterilization cycles with longer times than those listed are also acceptable.
- 8 Medical grade steam sterilization compatible wrap equivalent to four
- thicknesses of 140-thread-count muslin.
- Rigid sterilization container that compiles with ANSI/AAMI 5T46.
 Rash (unwrapped) sterilization by exposure at 132°C /270°F should only be used as an emergency procedure. Instruments must be cleaned and disassembled.
- 11 Drying times vary according to load size and should be increased for larger loads.
- 12 Gravity steam sterilization cycle parameters are available on request from Customer Service.

Note: The Sterillzer Manufacturer's instructions for operation and load configuration should be followed explicitly.



M. Storage Instructions

- Sterile, packaged instruments should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, and temperature/humidity extremes.
- Sterile instrument packages should be carefully examined prior to opening to ensure that package integrity has not been compromised.

Note: Maintenance of sterile package integrity is generally event related. If a sterile wrap is torn, perforated, shows any evidence of tampering or has been exposed to moisture, the instrument set must be repackaged and sterilized.

Note: If there is any evidence that the lid seal or filters on a sterilization container have been opened or compromised, the sterile filters must be replaced and the instrument set resterilized.

9. Hospital Responsibilities for Zimmer Loaner Sets

 Orthopaedic surgical instruments generally have a long service life; however, mishandling or inadequate protection can quickly diminish their life expectancy. Instruments which no longer perform properly because of long use, mishandling, or improper care should be returned to Zimmer to be discarded. Notify your Zimmer representative of any instrument problems.

- Loaner sets should undergo all steps of decontamination, cleaning, disinfection, inspection, and terminal sterilization before being returned to Zimmer. Documentation of decontamination should be provided with instruments being returned to Zimmer.
- Missing or damaged instruments from loaner sets should be brought to the attention of the operating room supervisor, to the director of the central supply department, and to your Zimmer representative to ensure that the next hospital will receive a complete set of instruments in working condition.
- The instructions provided in this manual have been validated by Zimmer in the laboratory and are capable of preparing orthopaedic devices for use. It is the responsibility of the Hospital to ensure that reprocessing is performed using the appropriate equipment and materials, and that personnel in the reprocessing facility have been adequately trained in order to achieve the desired result. Equipment and processes should be validated and routinely monitored. Any deviation by the processor from these instructions should be properly evaluated for effectiveness to avoid potential adverse consequences.

10. Customer Service Information

Mailing Address	Telephone				
Zimmer, Inc. 1800 West Center Street Warsaw, Indiana 46580 USA	Inside USA: 1-800-348-2759 Outside USA: local International access code +1-574-367-6131				
Zimmer GmbH Sulzer-Allee 8 CH-8404 Winterthur, Switzerland	+41 (0)52 262 60 70				
This Zimmer reprocessing manual and the associated Quick Reference Guide can be found at <u>www.zimmer.com</u> under the "Medical Professional" heading.					



11. References

- 1. AAMI TIR12, Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers
- 2. AAMI TIR1 3, Principles of industrial moist heat sterilization
- 3. AAMI TIR30. A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices
- ANSI/AAMI ST33, Guidelines for the selection and use of reusable rigid container systems for ethylene oxide sterilization and steam sterilization in health care facilities
- ANSI/AAMI ST35, Safe handling and biological decontamination of reusable medical devices in healthcare facilities and in nonclinical settings
- 6. ANSI/AAMI ST37, Flash sterilization Steam sterilization of patient care items for immediate use
- 7. ANSI/AAMI ST46, Steam sterilization and sterility assurance in health care facilities
- 8. ANSI/AAMI ST67, Sterilization of health care products Requirements for products labeled "Sterile"
- ANSI/AAMI STB1, Sterilization of medical devices Information to be provided by the manufacturer for the processing of resterilizable medical devices
- 10. ANSI/AAMI/ISO 15223 and Amendments 1 and 2, Medical devices Symbols to be used with medical device labels, labeling, and Information to be supplied
- 11. AORN, Standards, Recommended Practices and Guidelines
- 12. ASTM F 565, Standard Practice for Care and Handling of Orthopedic Implants and Instruments
- German Instrument Working Group (AKI) Arbeitskreis Instrumenten-Aufbereitung, Proper Maintenance of Instruments, 8th Ed, 2004.
- 14. IAHCSMM, Central Service Technical Manual
- 15. ISO 15883, Washer/Disinfectors: Requirements, Definitions and Test Methods
- ISO 17664, Sterilization of medical devices Information to be provided by the manufacturer for the processing of resterilizable medical devices

- Robert Koch Institute (RKI), Hospital Supplies and Instrument Sterilization in Light of CJD Patients and Suspected CJD Cases, Federal Health Gazette, 7/1998
- UK Department of Health, Health Technical Memorandum (HTM) 2010, Sterilization, Part 5 – Good Practice Guide
- 19. UK Department of Health, Health Technical Memorandum (HTM) 2030, Washer-Disinfectors – Validation and Verification
- 20. World Health Organization (WHO), WHO/CDS/CSR/APH 200.3, WHO Infection Control Guidelines for TSE

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017

Manual Orthopaedic Surgical Instruments

Appendix 1. Cleaning/Disinfection Processes

Chart 1. Manual Cleaning/Disinfection Procedure

Step 1	Completely submerge instruments in enzyme solution and allow to soak for 20 minutes. Scrub using a soft-bristled, mylon brush until all visible soil has been removed.
Step 2	Remove the device from the enzyme solution and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, holes and other difficult to reach areas.
Step 3	Place prepared cleaning agents in a sonication unit. Completely submerge device in cleaning solution and sonicate for 10 minutes at 45-50kHz.
Step 4	Rinse instrument in purified water for at least 3 minutes or until there is no sign of blood or soil on the device or in the rinse stream. Thoroughly and aggressively flush lumens, holes and other difficult to reach areas.
Step 5	Repeat the sonication and rinse steps above.
Step 6	Remove excess moisture from the instrument with a clean, absorbent and non- shedding wipe.

Chart 2. Combination Manual/Automated Cleaning/ Disinfection Procedure

Step 1	Completely submerge the Instruments in enzyme solution and allow to soak for 10 minutes. Use a soft, nyion-bristled brush to gently scrub the device until all visible soll has been removed. Particular attention must be given to crevices, lumens, mated surfaces, connectors and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft nyion-bristled brush.
Step 2	Remove devices from the enzyme solution and nose in purified water for a minimum of 1 minute. Thoroughly and aggressively flush tumens, holes and other difficult to reach areas.
Step 3	Place instruments in a suitable washer/disinfector basket and process through a standard washer/disinfector instrument cycle.

Chart 3. Typical Automated Washer/Disinfector Cycle for Cleaning/ Disinfection of Surgical Instruments

Step 1	Pre Wash; Cold Softened Tap Water; 2 minutes
Step 2	Enzyme Spray, Hot Softened Tap Water; 20 seconds
Step 3	Enzyme Soak; 1 minute
Step 4	Rinse (X2): Cold Softened Tap Water, 15 seconds
Step 5	Detergent Wash; Hot Softened Tap Water; (64-66°C/146-150°F); 2 minutes
Step 6	Rinse (X2); Hot Softened Tap Water: 15 seconds
Step 7	Thermal Rinse; Hot Softened Tap Water; (80-93°C/176-200°F); 2 minutes
Step 8	Punified Water Rinse; (64-66°C/146-150°F); 10 seconds
Step 9	Hot Air Dry; (116°C/240°F); 7 to 30 minutes

47



Manual Orthopaedic Surgical Instruments: Recommendations for Care, Cleaning, Maintenance and Sterilization

Page intentionally left blank.

Manual Orthopaedic Surgical Instruments: Recommendations for Care, Cleaning, Maintenance and Sterilization

Page intentionally left blank.

49

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017

Contact your Zimmer representative or visit us at www.zimmer.com



97-5000-170-00 Rev.3 2.5ML Printed in USA @1987, 2006, 2010 Zimmer, Inc. 50

3

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

Allen, Peter

From:	Hanley, Casey
Sent:	Friday, December 28, 2012 3:49 PM
To:	JOANNA.SURMA@ZIMMER:COM
Cc;	Allen, Peter
Subject:	K123862 - Refuse to Accept Notification
Attachments:	rtax_K123862.pdf

Dear Ms. Surma,

We have completed the administrative acceptance review of your premarket notification (510(k)) submission K123862) Our review indicates that your 510(k) submission does not meet the criteria established for administrative completeness.

Please refer to the attached checklist and e-mail your response referencing the 510(k) number K123862 by January 21, 2013. Your response should address all of the elements identified as missing or inconsistent in the attached checklist. Failure to provide the requested information by January 21, 2013 may result in placing the file on hold.

Upon receipt of the requested information, FDA may have additional requests for information.

Please be advised that once our RTA policy is implemented, if your 510(k) submission does not meet the criteria established for administrative completeness, your file will not be accepted and we will not begin our substantive review until you submit the missing elements and your submission is accepted. Refer to the draft <u>guidance</u> document for information.

Should you have questions about this email, you may either contact me or Peter Allen, the lead reviewer assigned to your 510(k) submission.

Sincerely, Casey

Casey L. Hanley, Ph.D. U.S. Food and Drug Administration CDRH/ODE/DOD/JFOB 10903 New Hampshire Ave WO66 Room 1567 Silver Spring, MD 20993 Tel: (301)796 6948 Casey.hanley@fda.hhs.gov

This e mail message is intended for the exclusive use of the recipient(s) named above. It may contain information that is protected, privileged, or confidential and it should not be disseminated distributed or copied to persons not authorized to receive such information. If you are not the intended recipient, any dissemination distribution or copying is strictly prohibited if you think you have received this e mail message in error please e-mail the sender immediately.

This communication is consistent with 21 CFR 10.85 (k) and constitutes an informal communication that represents my best judgment at this time but does not constitute an advisory opinion, does not necessarily represent the formal position of FDA, and does not bind or otherwise obligate or commit the agency to the views expressed.

Acceptance Checklist for Traditional 510(k)s

(should be completed within 15 days of DCC receipt)

The following information is not intended to serve as a comprehensive review.

510(k) Number: <u>K123862</u> Date Received: <u>December 14, 2012</u>

Lead Reviewer Name: <u>Casey Hanley</u> Branch: <u>JFOB</u> Division: <u>DOD</u> Office: <u>ODE</u>

Preliminary Questions				
Answers in the shaded blocks indicate consultation with Center advisor is needed.	Yes	No		
1. Is the product a device (per section 201(h) of the FD&C Act) or a combination product (per 21 CFR 3.2(e)) with a device constituent part subject to review in a				
510(k)?	x			
If it appears not to be a device (per section 201(h) of the FD&C Act) or such a combination product, or you are unsure, consult with the CDRH Jurisdictional Officer or the CBER Office Jurisdiction Liaison to determine the appropriate action, and inform division management. <i>Provide a summary of the Jurisdictional Officer's/Liaison's determination.</i> If the product does not appear to be a device or such a combination product, mark "No."				
Comments:				
2. Is the application with the appropriate Center? If the product is a device or a combination product with a device constituent part, is it subject to review by the Center in which the submission was received? If you believe the application is not with the appropriate Center or you are unsure, consult with the CDRH Jurisdictional Officer or CBER Office Jurisdiction Liaison to determine the appropriate action and inform your division management. <i>Provide a summary of the Jurisdictional</i> <i>Officer's/Liaison's determination</i> . If application should not be reviewed by your Center mark "No."	x			
Comments:				
3. Is a 510(k) the appropriate regulatory submission?				
If a 510(k) does not appear to be appropriate (e.g., Class III type and PMA required, or Class I or II type and 510(k)-exempt), you should consult with the CDRH 510(k) Program Director or appropriate CBER staff during the acceptance review. If 510(k) is not the appropriate regulatory submission, mark "No."	x			

Comments:	
 4. Is there a pending PMA for the same device with the same indications for use? If there is a pending PMA for the same device, consult division management and the CDRH 510(k) Program Director or appropriate CBER staff to determine the appropriate action. 	×
Comments:	
5. If clinical studies have been submitted, is the submitter the subject of an Application Integrity Policy (AIP)?	X
If yes, consult with the CDRH Office of Compliance/Division of Bioresearch Monitoring (OC/DBM - BIMO) or CBER Office of Compliance and Biologics Quality/Division of Inspections and Surveillance/Bioresearch Monitoring Branch (OCBQ/DIS/BMB) to determine the appropriate action. Check on web at http://www.fda.gov/ICECI/EnforcementActions/ApplicationIntegrityPolicy/ucm134453.ht	
<u>m</u> .	

If the answer to 1 or 2 appears to be "No," then stop review of the 510(k) and issue the "Original Jurisdictional Product" letter. If the answer to 3 is no, the lead reviewer should consult division management and other Center resources to determine the appropriate action.

If the answer to 4 is "Yes," then stop review of the 510(k), contact the CDRH 510(k) Staff and PMA Staff, or appropriate CBER staff.

If the answer to 5 is "Yes," then contact CDRH/OC/DBM – BIMO or CBER/OCBQ/DIS/BMB, provide a summary of the discussion with the BIMO Staff, and indicate BIMO's recommendation/action.

Organizational Elements

Failure to include these items alone generally should not result in an RTA designation

	Yes	No
a. Submission contains Table of Contents	X	
b. Each section is labeled (e.g., headings or tabs designating Device Description section, Labeling section, etc.)	x	
c. All pages of the submission are numbered All pages should be numbered in such a manner that information can be referenced by page number. This may be done either by consecutively numbering the entire submission, or numbering the pages within a section (e.g., 12-1, 12-2).	x	
d. Type of 510(k) is identified-traditional, abbreviated, or special If type of 510(k) is not designated, review as a traditional Comments:	x	

Elements of a Complete Submission (RTA Items) (21 CFR 807.87 unless otherwise indicated)

Submission should be designated RTA if not addressed

Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed. Any "No" answer will result in a "Refuse to Accept" decision. Yes N/A No Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission. Administrative A. 1. All content used to support the submission is written in English Х Π (including translations of test reports, literature articles, etc.) Comments: 2. 510(k) cover letter that identifies: Х a. Device trade name or proprietary name Х b. Device common name X П Х c. Device class and panel п Comments: Submission contains Indications for Use Statement with Rx and/or OTC 3. Х designated (see also 801.109) Submitter should use format appropriate for the reviewing Center/Office (CDRH/ODE, CDRH/OIVD, CBER/OBRR, CBER/OCTGT). If not provided in correct format, request the correct format during substantive review. Comments: 4. Submission contains 510(k) Summary or 510(k) Statement Х Π Either a) or b) must be answered "Yes" to be considered complete. Identify any missing element(s) as Comments. Summary contains all elements per 21 CFR 807.92 Ō $\overline{\Box}$ $\overline{\mathbf{X}}$ a. See also 510(k) Summary Checklist b. Statement contains all elements per 21 CFR 807.93 Comments: The 510(k) summary is missing a complete Device Description, (lacking

		<u>Elements of a Complete Submission (RTA Items)</u> (21 CFR 807.87 unless otherwise indicated)			
		Submission should be designated RTA if not addressed			
Check	"Yes	" if item is present, "N/A" if it is not needed and "No" if it is not includ	led bu	t neede	d.
-		 Any "No" answer will result in a "Refuse to Accept" decision. Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission. 	Yes	N/A	No
		explanation of how the device functions, the scientific concepts that f device, and the significant physical and performance characteristics device design, material used, and physical properties) and is missing predicate devices.	of the (device,	suc <u>h as</u>
	5.	Submission contains Truthful and Accuracy Statement per 21 CFR 807.87(k) See recommended format	х		
		Comments:			_
	6.	Submission contains Class III Summary and Certification See recommended <u>content</u> Form should be signed by a responsible person of the firm, not a consultant. CDRH is not currently able to accept a digital signature. "N/A" only if submission is not a Class III 510(k).		х	
		Comments:			
	7.	If submission relies upon a national or international standard as part of demonstration of substantial equivalence, submission contains Standards Data Report for 510(k)s (FDA Form 3654) or includes detailed information about how and the extent to which the standard has been followed There should be a completed form for each referenced national or international standard. "N/A" only if submission does not reference any standards.	x		
		Comments:			
	8.	Does submission contain clinical data? Select " N/A " for this item and 8.a. and 8.b. if the submission does not contain clinical data. If submission does contain clinical data, parts a. and b. must be answered "yes" for the $510(k)$ to be complete.		х	

			Elements of a Complete Submission (RTA Items) (21 CFR 807.87 unless otherwise indicated)						
		÷	Submission should be designated RTA if not addressed						
Checl	k "Yes	" if iten	n is present, "N/A" if it is not needed and "No" if it is not includ	led bu	t neede	:d.			
		 Any "No" answer will result in a "Refuse to Accept" decision. Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission. 							
		a.	Submission includes Financial Certification/Disclosure Statement		Х				
		b.	Submission includes Certification of Compliance with requirements of ClinicalTrials.gov Data Bank (FDA Form 3674) (42 U.S.C. 282(j)(5)(B))		х				
		Comm	nents:						
	9.	bundle device See <u>Gu</u> or Mu	is a bundled submission, the submission has been appropriately ed [i.e., the correct user fee(s) have been paid for the es/indications (section 738 of the FD&C Act)] uidance for Industry and FDA Staff: Bundling Multiple Devices htiple Indications in a Single Submission. ' if not a bundled submission		х				
205020		Comm	nents:		a .u				
	10,	which to sup Submi determ	Ibmission identifies prior submissions for the same device for FDA provided feedback related to the data or information needed port substantial equivalence (e.g., submission numbers for Pre- ission, IDE, prior not substantially equivalent (NSE) mination, prior 510(k) that was deleted or withdrawn) or states that were no prior submissions.	x					
		s r	f there were prior submissions: within current submission, the ponsor has identified where in the current submission any issues related to substantial equivalence outlined in prior communications are addressed.		x				
		Comm	nents:			_			
B	Dev	ice Des	cription						
	11.		e is a device-specific guidance document or special controls able to this submission, documentation has been provided to		x	۵			

.

<u>Elements of a Complete Submission (RTA Items)</u> (21 CFR 807.87 unless otherwise indicated)

Submission should be designated RTA if not addressed

heck '	"Yes"	" if iten	n is present, "N/A" if it is not needed and "No" if it is not inclu-	ded bu	t neede	đ
	•	Each subn any the c	"No" answer will result in a "Refuse to Accept" decision. In element on the checklist should be addressed within the mission. The submitter may provide a rationale for omission for criteria that are deemed not applicable. If a rationale is provided, criterion is considered present (Yes). An assessment of the male will be considered during the review of the submission.	Yes	N/A	No
		applica regard statuto Select omitte	ish that the submitter has followed the recommendations in the able device-specific guidance document or special controls ling the device description or otherwise met the applicable ory or regulatory criteria through an alternative approach. "No" if the submission does not include a rationale for any d information or any alternative approaches. "N/A" if there is no device-specific guidance document			
		Comm	ients:			
	12.	submi	scriptive information is present and consistent within the ssion (e.g., the device description section is consistent with the description in the labeling), including:			
	-		A description of the principle of operation and mechanism of action for achieving the intended therapeutic/diagnostic effect.	x		
	-	f	A description of all conditions of use such as surgical technique for implants; anatomical location of use; user interface; how the device interacts with other devices; and/or how the device interacts with the patient.	x		
		r	A list and description of each model for which clearance is requested. Select "N/A" if there is only one model.	x		
		Comm	nents:		_	
	13.	schem	e applicable, submission contains engineering drawing(s), atics, illustrations and/or figures of the device that are clear and e, and include dimensions	x		
			f engineering drawings, schematics, illustrations and/or figures are provided, one is provided for each model to be marketed	х		
		Comm	nents:	-		

			Elements of a Complete Submission (RTA Items) (21 CFR 807.87 unless otherwise indicated)							
	Submission should be designated RTA if not addressed									
Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.										
		Ar Ea sul an the rat	Yes	N/A	No					
	14.	acce Selei	evice is intended to be marketed with multiple components, ssories, and/or as part of a system, ct "N/A" if the device is not intended to be marketed with multiple ponents, accessories, and/or as part of a system.							
		a.	A description (as detailed in item 12.a. and b. and 13 above) is provided for each component or accessory.	x						
		b.	A 510(k) number is provided for each component or accessory that received a prior $510(k)$ clearance. Select "N/A" if the submission states that the component(s)/accessory(ies) do not have a prior $510(k)$ clearance.		х					
	1944 - 1944 - 1944 - 1944 - 1944 - 1944 - 1944 - 1944 - 1944 - 1944 - 1944 - 1944 - 1944 - 1944 - 1944 - 1944 -	Com	ments:							
C.	Sub	stanti	al Equivalence Discussion							
	15.	Subr	nitter has identified a predicate(s) device	x						
		a.	Predicate's 510(k) number, trade name, and model number (if applicable) provided	x						
		b.	The identified predicate(s) is consistent throughout the submission (i.e., the predicate(s) identified in the Substantial Equivalence section is the same as that listed in the 510(k) Summary (if applicable) and that used in comparative performance testing	x						
		Com	iments:							
	16.		nission includes a comparison of the following for the predicate(s) subject device							
		a.	Indications for use	X						
		b.	Technology, including features, materials, and principles of operation	х						

140

	Elements of a Complete Submission (RTA Items) (21 CFR 807.87 unless otherwise indicated)							
			Submission should be designated RTA if not addressed					
Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.								
-	 Any "No" answer will result in a "Refuse to Accept" decision. Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission. 				N/A	No		
		Com	iments:					
	17.	subj not c diffe <u>of th</u>	nission includes an analysis of why any differences between the ect device and predicate(s) do not render the device NSE (e.g., do constitute a new intended use, affect safety or effectiveness, or raise rent questions of safety and effectiveness) (see section 513(i)(1)(A) e FD&C Act) ere is no difference between the subject and predicate(s) with			X		
		resp state subn abov pote subn	ect to indications for use or technology, this should be explicitly ed, in which case "N/A" should be selected. Select "No" only if the nission does not include an analysis of differences as described be or a statement that there are no differences. Note that due to intial differences in manufacturing that may not be known to the nitter, no identified differences does not necessarily mean that no prmance testing is needed.					
		devi pred	ments: The sponsor provides a tabulated comparison of the subjects. From this tabulated list, it is clear there are differences betw licate devices with respect to indications for use and design feature not provide an analysis/discussion of why these differences do not provide an analysis/discussion of why these differences do not provide an analysis/discussion of why these differences do not provide an analysis/discussion of why these differences do not provide an analysis/discussion of why these differences do not provide an analysis/discussion of why these differences do not provide an analysis/discussion of why these differences do not provide an analysis/discussion of why these differences do not provide an analysis/discussion of why these differences do not provide an analysis/discussion of why these differences do not provide an analysis/discussion of why these differences do not provide an analysis/discussion of why these differences do not provide an analysis/discussion of why these differences do not provide an analysis/discussion of why these differences do not provide an analysis/discussion of why these differences do not provide an analysis/discussion of why these differences do not provide analysis/discussion diff	een th res. Th	e subje ie spon	ct a <u>nd</u> sor		
D.	Prop	osed	Labeling (see also 21 CFR part 801)					
	18.	pack	nission includes proposed labels, labeling (e.g., instructions for use, age insert, operator's manual), and advertisements that describe the ce, its intended use, and the directions for use	x				
		a.	Indications for use stated in labeling (21 CFR 801.61) and identical to IFU form and 510(k) Summary (if applicable)	х				
		b.	Directions for use included (including relevant hazards, warnings, precautions, contraindications), including directions for layperson (see 21 CFR 801.5) unless submission states that device qualifies	х		D		

٠,

			Elements of a Complete Submission (RTA Items) (21 CFR 807.87 unless otherwise indicated)						
Submission should be designated RTA if not addressed Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.									
	 Any "No" answer will result in a "Refuse to Accept" decision. Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission. 					No			
			for exemption per 21 CFR 801 Subpart D.		2				
		Com	ments:						
	19.	9. If indicated for prescription use, labeling includes the prescription use statement (see 21 CFR 801.109(b)(1)) or "Rx only" symbol [See also <u>Alternative to Certain Prescription Device Labeling Requirements</u>] Select "N/A" if not indicated for prescription use.							
		Com	Comments:						
	20.	Gene	ral labeling provisions						
		a.	Labeling includes name and place of business of the manufacturer, packer, or distributor (21 CFR 801.1)	x					
		b.	Labeling includes device common or usual name (21 CFR 801.61)	x					
		Com	ments:						
	21.	provi recor contr regul	If there is a device-specific guidance, special controls, or regulation, the provided labeling establishes that the submitter has followed the recommendations in the applicable guidance document, special controls, or regulation, or otherwise has met the applicable statutory or regulatory criteria through an alternative approach. Select "N/A" if there is no device-specific guidance or regulation.						
		Com	ments:						
	22.	all ap	e device is an in vitro diagnostic device, provided labeling includes oplicable information required per 21 CFR 809.10. et "N/A" if not an in vitro diagnostic device.		x				
Е.	Perf	ormai	nce Data – General						
	Submission: (one of the below must be checked) X does								

٩.

9

		Elements of a Complete Submission (RTA Items) (21 CFR 807.87 unless otherwise indicated)			
L		Submission should be designated RTA if not addressed			
Check "	Yes"	' if item is present, "N/A" if it is not needed and "No" if it is not includ	led bu	t neede	d.
	•		Yes	N/A	No
	conta	es not iin performance data. oes not" is selected, the performance data-related criteria below are omit	ted fro	m the ci	hecklist.
	Com	ments:			
	23.	Full test report is provided for each completed test to explain how the data generated from the test supports a finding of substantial equivalence. (A full test report includes: objective of the test, description of the test methods and procedures, study endpoint(s), pre-defined pass/fail criteria, results summary, and conclusions.)	x		
		Comments:			
	24.	Submission includes document to establish that the submitter has followed the recommendations for performance data outlined in the applicable device-specific guidance document or special controls, or otherwise met the applicable statutory or regulatory criteria through an alternative approach. Select "N/A" if there is no device-specific guidance document.	x		
		Comments:			
	25.	If literature was used as performance data, submission includes reprints or a summary of each article, and a discussion as to how each article is applicable to support the substantial equivalence of the subject device to the predicate.		x	
		Comments:			
	26.	If an animal study was conducted, Select "N/A" if no animal study was conducted.		х	
		a. Submission includes a study protocol and final study report to explain how the data generated from the study supports a finding of substantial equivalence. (A final study report includes a			

÷

Elements of a Complete Submission (RTA Items) (21 CFR 807.87 unless otherwise indicated)							
	_		Submission should be designated RTA if not addressed				
Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.							
	 Any "No" answer will result in a "Refuse to Accept" decision. Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission. 					No	
			contributing scientist report for each preclinical endpoint of evaluation within the protocol, such as performance/handling, in life, imaging, pathology, etc.)				
		ь.	Submission contains a statement that the study was conducted in compliance with applicable requirements in the GLP regulation (21 CFR Part 58), or, if the study was not conducted in compliance with the GLP regulation, the submission explains why the noncompliance would not impact the validity of the study data provided to support a substantial equivalence determination.				
		Com	ments:				
F.	Steri	ilizati	on				
	X sto X no I no	erile on-ster on-ster	n states that the device and/or accessories are: (<i>one of the below mus</i> ile but sterilized by the end user rile when used				
			nation will determine whether and what type of additional information for a substantial equivalence determination.	on may	v be		
	the c	heckli	erile when used" is selected, the sterility-related criteria below are o ist. tion regarding the sterility status of the device is not provided, select	-		·	
	Com	ments					
	27.	Asse	ssment of the need for sterilization information				
		a.	Identification of device, and/or accessories, and/or components that are provided sterile.	x			
		b.	Identification of device, and/or accessories, and/or components	ō	Ξ	$\overline{\mathbf{X}}$	

			<u>Elements of a Complete Submission (RTA Items)</u> (21 CFR 807.87 unless otherwise indicated)			
			Submission should be designated RTA if not addressed			
Check	"Yes	" if it	em is present, "N/A" if it is not needed and "No" if it is not include	led bu	t neede	d.
		Ea su an the	by "No" answer will result in a "Refuse to Accept" decision. Ich element on the checklist should be addressed within the bission. The submitter may provide a rationale for omission for y criteria that are deemed not applicable. If a rationale is provided, e criterion is considered present (Yes). An assessment of the tionale will be considered during the review of the submission.	Yes	N/A	No
			that are end user sterilized			
		с.	Identification of device, and/or accessories, and/or components that are reusable and cleaning/disinfection instructions are provided.			X
	Comments: The sponsor has identified a number of patient contacting instruments; however, the sponsor has not discussed whether they are provided sterile or non-sterile and whether they are able to be cleaned and re-sterilized. Therefore, please provide 1) a discussion as to which components are end user sterilized and 2) instructions for cleaning and re-sterilization for these components.					
		+				
	28.	Sele	e device, and/or accessory, and/or a component is provided sterile: ct "N/A" if no part of the device, accessories, or components is vided sterile, otherwise complete a-f below.		B	
	28.	Sele	e device, and/or accessory, and/or a component is provided sterile: ct "N/A" if no part of the device, accessories, or components is	x	B	
	28.	Sele prov	e device, and/or accessory, and/or a component is provided sterile: ct "N/A" if no part of the device, accessories, or components is ided sterile, otherwise complete a-f below. Sterilization method is stated for each component (including parameters such as dry time for steam sterilization, radiation dose,	x		
	28.	Selec prov a.	 e device, and/or accessory, and/or a component is provided sterile: ct "N/A" if no part of the device, accessories, or components is bided sterile, otherwise complete a-f below. Sterilization method is stated for each component (including parameters such as dry time for steam sterilization, radiation dose, etc.) A description of method to validate the sterilization cycle (e.g., half-cycle method and full citation of FDA-recognized standard, 			
	28.	Selec prov a. b.	 e device, and/or accessory, and/or a component is provided sterile: ct "N/A" if no part of the device, accessories, or components is nided sterile, otherwise complete a-f below. Sterilization method is stated for each component (including parameters such as dry time for steam sterilization, radiation dose, etc.) A description of method to validate the sterilization cycle (e.g., half-cycle method and full citation of FDA-recognized standard, including date) is provided For devices sterilized using chemical sterilants such as ethylene oxide (EO) and hydrogen peroxide, submission states maximum levels of sterilant residuals remaining on the device and sterilant residual limits. 	x		

.

			Elements of a Complete Submission (RTA Items) (21 CFR 807.87 unless otherwise indicated)			
			Submission should be designated RTA if not addressed			
Check	k "Yes	" if it	em is present, "N/A" if it is not needed and "No" if it is not includ	led bu	t neede	d
		Ea su an the	ny "No" answer will result in a "Refuse to Accept" decision. Ich element on the checklist should be addressed within the bmission. The submitter may provide a rationale for omission for y criteria that are deemed not applicable. If a rationale is provided, e criterion is considered present (Yes). An assessment of the tionale will be considered during the review of the submission.	Yes	N/A	No
		f.	If device is blood-contacting, a permanent implant, or contacts cerebrospinal fluid, or device is labeled "non-pyrogenic," submission contains a description of the endotoxin method used to make a determination (e.g., LAL), endotoxin release specification (e.g., 20 EU/device), and a rationale for the specification. Select "N/A" if device is not blood-contacting, not a permanent implant, does not contact cerebrospinal fluid, and is not labeled "non-pyrogenic." Select "N/A" if a rationale for omission is provided.			×
		prov	iments: (c) The sponsor has provided a promissory note for EO revided EO residuals remaining on the device following sterilization The sponsor does not provide a rationale for omission of Pyrogen	valida	tio <u>n.</u>	ie
	29.	guid indio regu	sterility information provided as recommended in the following ance documents or special controls, or information provided cating the submitter has otherwise met the applicable statutory or latory criteria through an alternative approach: er "a" or "b" must be answered "Yes" to be considered complete.			
		a.	Device-specific guidance document or special controls Select "N/A" if no device-specific guidance document.		х	
		b.	Cross-cutting guidance document (for more information see "Submission and Review of Sterility Information in Premarket Notification (510(k)) Submissions for Devices Labeled as Sterile") Select "N/A" if device-specific guidance followed instead.	X		D
		Соп	iments:			
G.	She	f Life				
	30.		e device is provided sterile or the device is provided non-sterile and age conditions (i.e., aging) could impact device safety or			

ī,

Contains Nonbinding Recommendations Draft - Not for Implementation

			<u>Elements of a Complete Submission (RTA Items)</u> (21 CFR 807.87 unless otherwise indicated)			_				
			Submission should be designated RTA if not addressed							
Chec	k "Yes	" if it	em is present, "N/A" if it is not needed and "No" if it is not inclu	ded bu	t neede	d				
		Ea su an the	by "No" answer will result in a "Refuse to Accept" decision. Inch element on the checklist should be addressed within the bission. The submitter may provide a rationale for omission for y criteria that are deemed not applicable. If a rationale is provided, e criterion is considered present (Yes). An assessment of the ionale will be considered during the review of the submission.	Yes	N/A	No				
		Sele state	ctiveness, address the following: ct "N/A" if the device is not provided sterile and the submitter cs that storage conditions could not affect device safety or ctiveness.							
		a.	Proposed shelf life/expiry date stated	x						
		b.	Submission includes description of shelf life validation method(s) used to ensure device performance and sterility, as applicable, remain substantially equivalent to that of the predicate device throughout the stated shelf life.	x						
		Com	iments:	.e						
H.	Bioc	ompa	atibility							
	X and ardirect This nece	Submission states that there: (<i>one of the below must be checked</i>) X are are not direct or indirect (e.g., through fluid infusion) patient-contacting components. This information will determine whether and what type of additional information may be necessary for a substantial equivalence determination.								
	chec selec	klist. ct "No				·				
	Com	ment:	······································							
	31.	asso	nission includes list of patient-contacting device components and ciated materials of construction, including identification of color tives, if present	X						
ur denduly f		Com	iments:							

.

.

.

Elements of a Complete Submission (RTA Items) (21 CFR 807.87 unless otherwise indicated)								
Submission should be designated RTA if not addressed								
Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.								
	 Any "No" answer will result in a "Refuse to Accept" decision. Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission. 							
	32. Submission identifies contact classification (e.g., surface-contacting, less than 24 hour duration)							
_		Comments: Sponsor does not state the contact classification for the in	strum	ents.				
	33.	Biocompatibility assessment of patient-contacting components Submission includes: Test protocol (including identification and description of test article), methods, pass/fail criteria, and results provided for each completed test,	x					
		OR a statement that biocompatibility testing is not needed with a rationale (e.g., materials and manufacturing/processing are identical to the predicate).						
Ι.	Soft	ware			_			
	□ do X do	nission states that the device: (<i>one of the below must be checked</i>) bes bes not ain software.						
		information will determine whether and what type of additional informations sary for a substantial equivalence determination.	on may	/ be	P			
		oes not" is selected, the software-related criterion is omitted from the cher mation regarding whether the device contains software is not provided, se						
	Com	ments:						
	34.	All appropriate categories of software verification and validation documentation provided based on stated level of concern, as described in <u>Guidance for the Content of Premarket Submissions for Software</u> <u>Contained in Medical Devices</u> , or the submitter has provided						

.

<u>Elements of a Complete Submission (RTA Items)</u> (21 CFR 807.87 unless otherwise indicated)								
		Submission should be designated RTA if not addressed						
Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.								
-		 Any "No" answer will result in a "Refuse to Accept" decision. Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission. 						
		documentation that it has otherwise met the applicable statutory or regulatory criteria through an alternative approach.						
		Comments:						
J.	EMC and Electrical Safety							
	Submission states that the device: (one of the below must be checked) does X does not require EMC and Electrical Safety evaluation.							
	This information will determine whether and what type of additional information may be necessary for a substantial equivalence determination. If "does not" is selected, the EMC-related and Electrical Safety-related criteria below are							
	omit	ted from the checklist. If information regarding whether the device require trical Safety evaluation is not provided, select "No."						
	Com	ments:						
	35.	Submission includes evaluation of electrical safety per IEC 60601-1 or equivalent FDA-recognized standard and if applicable, the device- specific standard OR submission includes electrical safety evaluation using methods or standards that are not FDA-recognized and information indicating that these methods/standards otherwise meet applicable statutory and regulatory requirements.						
		Comments:						
	36.	Submission includes evaluation of electromagnetic compatibility per IEC 60601-1-2 or equivalent FDA-recognized standard and if applicable, the device-specific standard						

			<u>Elements of a Complete Submission (RTA Items)</u> (21 CFR 807.87 unless otherwise indicated)					
Submission should be designated RTA if not addressed								
Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.								
	•	Ea sui an the	by "No" answer will result in a "Refuse to Accept" decision. Inch element on the checklist should be addressed within the bission. The submitter may provide a rationale for omission for y criteria that are deemed not applicable. If a rationale is provided, e criterion is considered present (Yes). An assessment of the tionale will be considered during the review of the submission.	Yes	N/A	No		
		meth indic	nission includes electromagnetic compatibility evaluation using nods or standards that are not FDA-recognized and information cating that these methods/standards otherwise meet applicable atory and regulatory requirements.					
		Com	iments:					
К.			nce Characteristics – In Vitro Diagnostic Devices Only (see also 09.10(b)(12))					
	Submission indicates that device: (one of the below must be checked) is X is not an in vitro diagnostic device (IVD). If "is not" is selected, the performance data-related criteria below are omitted from the checklist.							
		ment						
	37.		nission includes the following analytical studies, including ciated protocols and line data:					
		a.	Precision/reproducibility (at least 3 sites generally necessary)					
		b.	Accuracy (includes linearity, assay cut-off, method comparison or comparison to clinical outcome, matrix comparison, reference range, and stability protocol and acceptance criteria)					
		c.	Sensitivity (detection limits (LoB, LoD, and LoQ))					
		d.	Analytical specificity					

<u>Elements of a Complete Submission (RTA Items)</u> (21 CFR 807.87 unless otherwise indicated)

Submission should be designated RTA if not addressed

Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.

 Any "No" answer will result in a "Refuse to Accept" decision. Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission. 	Yes	N/A	No
 Comments:			

Decision: Accept ____ Refuse to Accept _X___

3

If Accept, notify applicant; if Refuse to Accept, notify applicant in writing and include a copy of this checklist.

Digital Signature Concurrence Table			
Reviewer Sign-Off	Casey Hanley 2012.12.27_09:41:57 -05'00'		
Branch Chief Sign-Off	Anton E. Dmitriev 2012.12.27 14:45:32 -05'00'		
Division Sign-Off	Mark N. Melkerson 2012.12.27 15:30:01 -05'00'		

19



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

U.S. Food and Drug Administration Center for Devices and Radiological Health Document Control Center WO66-G609 10903 New Hampshire Avenue Silver Spring, MD 20993-0002

December 14, 2012

ZIMMER, INC. 1800 WEST CENTER STREET WARSAW, INDIANA 46580 ATTN: JOANNA L, SURMA 510k Number: K123862 Received: 12/14/2012 Product: ZIMMER NEXEL TOTAL ELBOW

The Food and Drug Administration (FDA), Center for Devices and Radiological Health (CDRH), has received the Premarket Notification, (510(k)), you submitted in accordance with Section 510(k) of the Federal Food, Drug, and Cosmetic Act(Act) for the above referenced product and for the above referenced 510(k) submitter. Please note, if the 510(k) submitter is incorrect, please notify the 510(k) Staff immediately. We have assigned your submission a unique 510(k) number that is cited above. Please refer prominently to this 510(k) number in all future correspondence that relates to this submission. We will notify you when the processing of your 510(k) has been completed or if any additional information is required. YOU MAY NOT PLACE THIS DEVICE INTO COMMERCIAL DISTRIBUTION UNTIL YOU RECEIVE A LETTER FROM FDA ALLOWING YOU TO DO SO.

Please remember that all correspondence concerning your submission MUST be sent to the Document Mail Center (DMC) at the above letterhead address. Correspondence sent to any address other than the one above will not be considered as part of your official 510(k) submission.

On September 27, 2007, the President signed an act reauthorizing medical device user fees for fiscal years 2008 - 2012. The legislation - the Medical Device User Fee Amendments of 2007 is part of a larger bill, the Food and Drug Amendments Act of 2007. Please visit our website at

http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/Overview/MedicalDeviceUserFeeandMod ernizationActMDUFMA/default.htm

for more information regarding fees and FDA review goals. In addition, effective January 2, 2008, any firm that chooses to use a standard in the review of ANY new 510(k) needs to fill out the new standards form (Form 3654) and submit it with their 510(k). The form may be found at <u>http://www.fda.gov/AboutFDA/ReportsManualsForms/Forms/default.htm</u>.

We remind you that Title VIII of the Food and Drug Administration Amendments Act of 2007 (FDAAA) amended the PHS Act by adding new section 402(j) (42 U.S.C. § 282(j)), which expanded the current database known as ClinicalTrials.gov to include mandatory registration and reporting of results for applicable clinical trials of human drugs (including biological products) and devices. Section 402(j) requires that a certification form <u>http://www.fda.gov/AboutFDA/ReportsManualsForms/Forms/default.htm</u> accompany 510(k)/HDE/PMA submissions. The agency has issued a draft guidance titled: "Certifications To Accompany Drug, Biological

Product, and Device Applications/Submissions: Compliance with Section 402(j) of The Public Health Service Act, Added By Title VIII of The Food and Drug Administration Amendments Act of 2007" <u>http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/HowtoMarketYourDevice/PremarketSubmissions/PremarketNotification510k/ucm134034.htm</u>. According to the draft guidance, 510(k) submissions that do not contain clinical data do not need the certification form.

Please note the following documents as they relate to 510(k) review: 1) Guidance for Industry and FDA Staff entitled, "Interactive Review for Medical Device Submissions: 510(k)s, Original PMAs, PMA Supplements, Original BLAs and BLA Supplements". This guidance can be found at

<u>http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm089402.htm</u>. Please refer to this guidance for information on a formalized interactive review process. 2) Guidance for Industry and FDA Staff entitled, "Format for Traditional and Abbreviated 510(k)s". This guidance can be found at <u>http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ ucm084365.htm</u>. Please refer to this guidance for assistance on how to format an original submission for a Traditional or Abbreviated 510(k).

In all future premarket submissions, we encourage you to provide an electronic copy of your submission. By doing so, you will save FDA resources and may help reviewers navigate through longer documents more easily. Under CDRH's e-Copy Program, you may replace one paper copy of any premarket submission (e.g., 510(k), IDE, PMA, HDE) with an electronic copy. For more information about the program, including the formatting requirements, please visit our web site at

http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/HowtoMarketYourDevice/PremarketSubmissio ns/ucm134508.html. In addition, the 510(k) Program Video is now available for viewing on line at http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/HowtoMarketYourDevice/PremarketSubmissio ns/PremarketNotification510k/ucm070201.htm.

Please ensure that whether you submit a 510(k) Summary as per 21 CFR 807.92, or a 510(k) Statement as per 21 CFR 807.93, it meets the content and format regulatory requirements.

Lastly, you should be familiar with the regulatory requirements for medical devices available at Device Advice <u>http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/default.htm</u>. If you have questions on the status of your submission, please contact DSMICA at (301)796-7100 or the toll-free number (800)638-2041, or at 'heir internet address <u>http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/default.htm</u>. If you have procedural questions, please contact the 510(k) Staff at (301)796-5640.

Sincerely,

510(k) Staff

Pugh, Dominique *

oft Outlook
.Surma@Zimmer.com
December 14, 2012 4:35 PM
I: K123862 ACK LETTER & ECOPY ATTACHMENT

Delivery to these recipients or groups is complete, but no delivery notification was sent by the destination server:

Joanna.Surma@Zimmer.com (Joanna.Surma@Zimmer.com)

Subject: K123862 ACK LETTER & ECOPY ATTACHMENT

EXPLANATION FOR ECOPY ATTACHMENT

You provided an eCopy along with the required number of paper copies. Your eCopy did not pass our validation software. The specific reasons for this failure are identified in the attachment.

Because eCopies are currently voluntary and because you already provided the necessary number of paper copies, you do <u>NOT</u> need to provide a replacement eCopy to FDA.

Instead, FDA is providing you with this information so that you can understand the reasons why your eCopy failed in order to better assure that an eCopy that you submit for another submission will not have similar issues.

As additional resources:

- Refer to our eCopy guidance at <u>http://www.fda.gov/downloads/MedicalDevices/DeviceRegulationand</u> <u>Guidance/GuidanceDocuments/UCM313794.pdf</u>. The eCopy program is not implemented (i.e., eCopy remain voluntary) at this time; however, the technical specifications/standards in Attachment 1 of the guidance are already in place.
- We strongly encourage you to take advantage of our free new eSubmitter-eCopies tool at <u>http://www.fda.gov/ForIndustry/FDAeSubmitter/ucm317334.htm</u>, which will help create an eCopy that passes the new technical specifications/standards.

Please keep in mind that once we implement the new eCopy program, most submissions types will go on eCopy hold until a valid eCopy is received.

If you have any questions about the eCopy program or the technical specifications/standards, please contact Ms. Samie Allen at <u>samie.allen@fda.hhs.gov</u> or at 301-796-6055.

Attachment

Attachment for Sub	mission Numb	er(s):				
K123862		<u>ः</u>		and a second s	ing internet Series and series and s	2 1 1 1 1 1
		18 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	an than an a			i king i

The list below identifies the reason(s) why your eCopy failed FDA's eCopy validation process. All of these items need to be addressed or your eCopy will not pass the validation process.

1. The following PDF file(s) exceeded maximum allowed file size limit of 50MB:

002_C-20121201_2012-12-13_001_02_FDA_Submission.pdf

¢

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017



1.

V 1 St FDA CDRH DMC DEC 1 4 2012 Received

P.O. Box 708 Warsaw, IN 46581-0708 (574) 267-6131

13 December 2012

U.S. Food and Drug Administration Center for Devices and Radiological Health Document Mail Center – W066-G609 10903 New Hampshire Avenue Silver Spring, MD 20993-0002

Dear Sir or Madam:

Subject: Traditional 510(k) Premarket Notification – Zimmer Nexel Total Elbow – Prosthesis, Elbow, Constrained, Cemented – JDC – Class II – 21 CFR 888.3150

As required by Section 510(k) of the Federal Food, Drug and Cosmetic Act, as amended by the Medical Device Amendments of 1976, the Safe Medical Devices Act of 1990, FDA Modernization Act of 1997 (FDAMA) and in accordance with Title 21 of the Code of Federal Regulations (CFR) Part 807, subpart E, the above noted premarket notification is hereby submitted to the Food and Drug Administration (FDA). As required by 21 CFR 807.90(c), this document is submitted in duplicate. A CD-ROM containing an eCopy that is an exact duplicate of the paper copy submission is also included for your convenience.

These devices have not been previously submitted to FDA for identical or different indications, are not currently being reviewed for different indications by the same or different branch within ODE, and have not been previously cleared by FDA for different indications.

If you require any additional information or have any questions, please contact me by telephone at (574) 371-1642, e-mail at Joanna.Surma@Zimmer.com or fax at (574) 372-4605.

Sincerely,

form I durm 12 Dec 2012

Joanna L. Surma Associate Project Manager, Regulatory Affairs

Enclosure

٩.

14

÷

ŧ

	orni Approved, OMB No. 0910-5111 spiratum Date: February 28, 2013. See Instructions for OMB Statement			
DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION MEDICAL DEVICE USER FEE COVER SHEET	PAYMENT IDENTIFICATION NUMBER: (b)(4) Write the Payment Identification number on your check.			
A completed cover sheat must accompany each original application or supplement subject to fees. If payment is sent by U.S. mail or courier, please include a copy of this completed form with payment. Payment and mailing instructions can be found at http://www.fda.gov/oc/mdufma/coversheet.html				
1. COMPANY NAME AND ADDRESS (include name, street address, city state, country, and post office code)	2. CONTACT NAME Joanna Surma			
	2.1 E-MAIL ADDRESS			
	Joanna.Surma@Zimmer.com			
345 EAST MAIN STREET WARSAW IN 46580	2.2 TELEPHONE NUMBER (include Area code)			
US	574-371-1642			
1.1 EMPLOYER IDENTIFICATION NUMBER (EIN) *****5416	2.3 FACSIMILE (FAX) NUMBER (Include Area code)			
3. TYPE OF PREMARKET APPLICATION (Select one of the follow descriptions at the following web site: http://www.fde.gov/oc/mdufma	ing in each column; if you are unsure, please refer to the application			
Select an application type:	3.1 Select a center			
[X] Premarket notification(510(k)); except for third party	(X) CDRH			
[] 513(g) Request for Information	[]CBER			
[] Biologics License Application (BLA)	3.2 Select one of the types below			
[] Premarket Approval Application (PMA)	[X] Original Application			
[] Modular PMA	Supplement Types:			
[] Product Development Protocol (PDP)	[] Efficacy (BLA)			
[] Premarket Report (PMR)	[] Panel Track (PMA, PMR, PDP)			
[] Annual Fee for Periodic Reporting (APR)	() Real-Time (PMA, PMR, PDP)			
[] 30-Day Notice	[] 180-day (PMA, PMR, PDP)			
 4. ARE YOU A SMALL BUSINESS? (See the instructions for more information on determining this status) [] YES, I meet the small business criteria and have submitted the required [X] NO, I am not a small business qualifying documents to FDA 4.1 If Yes, please enter your Small Business Decision Number: 				
5. FDA WILL NOT ACCEPT YOUR SUBMISSION IF YOUR COMP THAT IS DUE TO FDA. HAS YOUR COMPANY PAID ALL ESTABL (X) YES (All of our establishments have registered and paid the fee	ISHMENT REGISTRATION FEES THAT ARE DUE TO FDA?			
 [X] YES (All of our establishments have registered and paid the fee, or this is our first device, and we will register and pay the fee within 30 days of FDA's approval/clearance of this device) [] NO (If "NO," FDA will not accept your submission until you have paid all fees due to FDA. This submission will not be processed; see 				
http://www.fda.gov/cdr/mdufma for additional information)				
6. IS THIS PREMARKET APPLICATION COVERED BY ANY OF T APPLICABLE EXCEPTION.	~			
[] This application is the first PMA submitted by a qualified small business. [] The sole purpose of the application is to support including any affiliates conditions of use for a pediatric population				
[] This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only				
7. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION FOR WHICH FEES WERE WAIVED DUE TO SOLE USE IN A PEDIATRIC POPULATION THAT NOW PROPOSES CONDITION OF USE FOR ANY ADULT POPULATION? (If so, the application is subject to the fee that applies for an original premarket approval application (PMA), [] YES [X] NO				
PAPERWORK REDUCTION ACT STATEMENT Public reporting burden for this collection of information is estimated to average 18 minutes per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to the address below.				
Department of Health and Human Services, Food and Drug Administration, Office of Chief Information Officer, 1350 Piccard Drive, 4th Floor Rockville, MD 20850 Please do NOT return this form to the above address, except as it pertains to comments on the burden estimate.]				
8. USER FEE PAYMENT AMOUNT SUBMITTED FOR THIS PREM	ARKET APPLICATION			
	30-Nov-2012			

"Close Window" Print Cover sheet

	DEPARTMENT OF HEALTH AND FOOD AND DRUG ADM		ES		Form Appro OMB No. 0	910-0120			
CDRH PREMARKET REVIEW SUBMISSION COVER SHEET Expiration Date: December 31, 2013 See OMB Statement on page 5.									
Date of Submission	User Fee Payment	ID Number		FDA Submiss	ion Docume	nt Numbe	er (if known)		
12/13/2012	(b)(4)								
SECTION A		TYPE OF SI	UBMISSION						
PMA	PMA & HDE Supplement	PD	P	510(K)	.		Meeting		
Original Submission	Regular (180 day)	Original PC		Criginal Subm	ission [.]	=	-510(K) Meeting		
Strengtheright Premarket Report Modular Submission	Special	Notice of C		X Traditional		·	-IDE Meeting PMA Meeting		
Amendment	30-day Supplement	Arnendmer		Special Abbreviated	(Complete		-PDP Meeting		
Report	30-day Notice	1		Abbreviated section I, Pa	age 5)	<u> </u>	y 100 Meeting		
Report Amendment	135-day Supplement			Additional Info	rmation	🗍 Agi	reement Meeting		
Licensing Agreement	Real-time Review			Third Party			termination Meeting		
	HDE Supplement	Į				01	er (specify):		
	Other	ł	l						
IDE	Humanitarian Device Exemption (HDE)	Class II Exemp	otion Petition	Evaluation of A Class III Desig	nation	Oth	er Submission		
[] Onginal Submission	Original Submission	🗌 Öriginal Su	bmission	(De Novo		513)(g)		
	Amendment	Additional	information	Additional Info	1	Oth (de	er scribe submission}:		
Supplement	Supplement	ļ			[100	achie aubmission;		
	Report Report Amendment								
Have you used or cited Stan		_ Yes [_]N⊄		please complete Se	ection I, Pagi	e 5)			
SECTION B Company / Institution Name		ITTER, APPLI	Establishment	Registration Number	(if known)				
Zimmer, Inc.									
Division Name (if applicable)			Phone Number	(including area code)				
			(574) 371 - 16-	12					
Street Address			FAX Number (i	ncluding area code)					
1800 West Center Street			(574) 372 - 460	er					
City			State / Province		ZIP/Postal	Code	Country		
Warsaw			Indiana	-	46580		USA		
Contact Name		·				·	L		
Joanna L. Surma									
Contact Title			Contact E-mail	Address					
Associate Project Manager, Re	gulatory Affairs			@Zimmer.com					
SECTION C		BONDENT /			n about				
Company / Institution Name	APPLICATION CORRES	I ONDENT (e.	g., consultan	, a onterent mor	n auove)				
Division Name (if applicable)			Phone Number	(including area code	}				
Street Address			FAX Number (ii	ncluding area code)			<u> </u>		
			······························	•					
City			State / Province		ZIP Code		Country		
Contact Name		· · · · · · · · · · · · · · · · · · ·	. <u> </u>		<u> </u>		L		
			A						
Contact Title			Contact E-mail	Address					
			L		_		·····		
FORM FDA 3514 (12/10)							age 1 of 5 Pages		

a,

2

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118 2

SECTION D1 REA	SON FOR APPLICATION - PMA, PDP, OR I	105
New Device Withdrawal Additional or Expanded Indications Request for Extension Post-approval Study Protocol Request for Applicant Hold Request for Removal of Applicant Hold Request to Remove or Add Manufacturing Site Process change: Manufacturing Packaging Sterilization Other (specify below) Response to FDA correspondence:	Son Fox Application - PMA, PDP, OR Change in design, component, or specification Software / Hardware Color Additive Material Specifications Other (specify below) Labeling change: Indications Instructions Performance Charactenstics Shelf Life Trade Name Other (specify below)	Location change. Manufacturer Sterilizer Packager Report Submission: Annual or Periodic Post-approval Study Adverse Reaction Device Defect Amendment Change in Ownership Change in Correspondent Change of Applicant Address
Other Reason (specify):		
SECTION D2 New Device New Indication Addition of Institution Expansion / Extension of Study IRB Certification Termination of Study Withdrawal of Application Unanticipated Adverse Effect Notification of Emergency Use Compassionate Use Request Treatment IDE Continued Access	REASON FOR APPLICATION - IDE Change in: Correspondent/Applicant Design / Device Informed Consent Manufacturer Manufacturing Process Protocol - Feasibility Protocol - Other Sponsor Current Investigator Annual Progress Report Site Waiver Report Final	 Response to FDA Letter Concerning Conditional Approval Deerned Approved Deficient Final Report Deficient Progress Report Deficient Investigator Report Disapproval Request Extension of Time to Respond to FDA Request Hearing
Other Reason (specify):	REASON FOR SUBMISSION - 510(k)	
X New Device	Additional or Expanded Indications	Change in Technology
Other Reason (specify):		Page 2 of 5 Pages

Ŧ

i

h

SE	CTIONE			ADDIT	ю	NAL INFORMATIO	N'ON 51	0(K	() SUB	MIS	sio	NS			
Pro	duct codes of device	s to v	hich substantia	l equivalen)Ce	IS claimed		-						or statement con ectiveness inform	
1	JDC		2			3	2	1						k) summary atta	
5			6			7	1	i i						k) statement	
Info	ormation on devices to	o whi	ch substantial e	auvalence	e i e	claimed <i>(if known</i>)	N		_				<u> </u>		
,			lumber	<u> </u>	2	Trade or Proprie	etary or M	odel	Name		2.		Ма	nulacturer	
	K001918		_			the Coonrad/Morrey To	otal Elbow	, ,	-			Zin	mer		
1]	1	-					11				
					_		· · · <u></u>				╢╌┥				
2	K053189				2	the Coonrad/Morrey T	otal Elboy	12 C			2	Zim	uner		
						· · · · · · · · · · · · · · · · · · ·									
3				l)	3						3				
3					3						3				
						· · · · · · · · · · · · · · · · · · ·									
4					4						4				
╞─					_						+			·	
5				- II	5						5				
									_						
6					6						6				
					-					_	1				
	mmon of usual name osthesis. Elbow, Cons			1e											
			- <u></u>							- 1					
	Trade or Proprietary	/ or M	odel Name for	This Device	ę	<u></u>				M	del	lint	er		
1	Zimmer Nexel Tota	1 Elbo	w						- II -	1 N	A				
2				• • •		<u> </u>				2					
Ľ.		• • •			_		·		'	·					
з									:	3					
										Ť					
4			_							4				<u>.</u>	2 <u>-</u> 2
5										5					
FD	A document numbers	s of al	l prior related s	ubmissions	s (re	gardless of outcome)							<u> </u>		
1			2		3		4			-	5		· · · ·	6	
7			8		9		10		-		11			12	
	ta Included in Submis								_				<u> </u>		
Da	a maadeo in Suomis	a si Ul I	🔀 Lai	boratory Te	stir	ng 📋	Animal Tri	als					Human Trials		
	CTION G				AS	SIFICATION - APP	LICATI	ON				ICA'	TIONS		
	duct Code		R. Section (if a	pplicable)					Device	Clas	s				
	XC	21	CFR 888.3150			· · · · · · · · · · · · · · · · · · ·				Class	I	×	Class II		
	ssification Panel									Class	ш	E	Unclassified		
0	rthopedic														
In cc or (6	ntributing to elbow in degenerative arthritis	ide: (1 istabil with articul	ity, (3) Ancient incapacitating p lar fracture of th	sepsis that min, (5) Ins we elbow joi	res stab	ch significantly comprov ults in ankylosed joints, ility or loss of motion w urfaces that precludes he ted use only.	especially hen the de	in ca gree	ases of b of joint	oilater or so	al an ft tiss	kylos ue da	is, (4) Advanced r mage precludes re	heumatoid, post- hable osteosynti	traumatic. resis,
	RM FDA 3514 (1)	2/101		···-							-			Page 3 of 5	Pages

			FDA Document Number (if kn	iown)	
Note: Submission of the inneed to submit device est	information entered in Section H do tablishment registration.	es not affect the			
SECTION H			TERILIZATION SITES RE	LATING TO A S	UBMISSION
X Original	Facility Establishment Identifier (I	FEI) Number	X Manufacturer	Contract Ster	alizer
Add Delete			Contract Manufacturer	🔲 Repackager	/ Relabeler
Company / Institution Nat	me		Establishment Registration N	umber	
Zimmer, Inc			1822565		
Division Name (if applicat	ble)	0-11	Phone Number (including are	a code)	
			(574) 372-4964		
Street Address			FAX Number (including area	code)	
1800 West Center Street	ι		(574) 372-4605		
City	· · · ·		State / Province	ZIP Code	Country
Warsaw			Indiana	46580	USA
Contact Name		Contact Title		Contact E	-mail Address
Carol Vierling		Director, Corporate	Regulatory Affairs	Carol Vi	ierling@Zimmer.com
	Facility Establishment Identifier (FEI) Number	Manufacturer	X Contract Ste	ohaor
Original			Contract Manufacturer	Repackager	
Company / Institution Nar					
(b)(4)			Establishment Registration N	unidei	
Drvision Name (if application	ble)	·······	Phone Number (Including are	a code)	
			(b)(4)		
Street Address	·		FAX Number (including area	code)	
(b)(4)					
City			State / Province	ZIP Code	Country
(b)(4)					
Contact Name	<u> </u>	Contact Title		Contact É	-mail Address
(b)(4)		oonaet miç		Combot E	
X Original	Facility Establishment Identifier (FEI) Number	Manufacturer	X Contract Ste	nlizer
Add Delete			Contract Manufacturer	Repackager	
Company / Institution Nar	 me		Establishment Registration N	umber	
(b)(4)					
Division Name (if applica	ble)		Phone Number (including are	a code)	
(b)(4)					
Street Address			FAX Number (including area	code)	
(b)(4)					
City			State / Province	ZIP Code	Country
(b)(4)					
Contact Name		Contact Title		Contact E	-mail Address
(b)(4)					

FORM FDA 3514 (12/10)

i

Add Continuation Page: Page 4 of 5 Pages

SE	СT	10	N	1
	-			

	Olevada Ma		Chandrade Tale	The set of the set	
	Standards No. 14283	Standards Organization ISO	Standards Title Implants for surgery – Fundamental principles	Version	Date
1					07/08/2004
	Standards No.	Standards Organization ISO	Standards Title Non-active surgical implants General requirements	Version	Date
2	24				01/07/2008
	Standards No. 11980	Standards Organization ASTM	Standards Title Standard Guidance for Accelerated Aging of Sterile Barrier Systems for Medical Devices	Version	Date
3					04/01/2007
	Standards No. ST81	Standards Organization	Standards Title Sterilization of medical devices - Information to be provided by the manufacturer for the processing of resterilizable medical devices	Version	Date
4			manufacturer for the processing of resternizable medical devices		12/02/2004
	Standards No. 10993-1	Standards Organization ANSI	Standards Title Biological Evaluation of Medical Devices – Part 1: Evaluation and	Version	Date
5		AAMI ISO	Testing within a risk management process		10/15/2009
	Standards No. 10993-3	Standards Organization	Standards Title Biological Evaluation of Medical Devices – Part 3: Tests for	Version	Date
6		ANSI AAMI ISO	genoloxicity, carcinogenicity and reproductive loxicity		10/15/2003
	Standards No.	Standards Organization ANSI	Standards Title Biological Evaluation of Medical Devices – Part 5: Tests for in vitro	Version	Date
7		AAMI ISO	cytotoxicity		06/01/2009
		Pleas	e include any additional standards to be cited on a separate pag	l	[

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer 1350 Piccard Drive, Room 400 Rockville, MD 20850

An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMR control number

FORM FDA 3514 (12/10)

Page 5 of 5 Pages

SECT	FION	1
------	-------------	---

UTILIZATION OF STANDARDS

	Standards No. 10993-6	Standards Organization ANSI AAMI ISO	Standards Title Biological Evaluation of Medical Devices – Part 5: Tests for in vitro cytoroxicity	Version	Date 4/15/2007
2	Standards No. 10993-7	Standards Organization ANSI AAMI ISO	Standards Title Biological Evaluation of Medical Devices – Part 7. Ethylene oxide sterilization residuals	Version	Date 10/15/2008
3	Standards No. 10993-10	Standards Organization ANSI AAMI ISO	Standards Title Biological Evaluation of Medical Devices – Part 10, Tests for irritation and delayed-type hypersensitivity Amendment 1	Version	Date 8/01/2006
4	Standards No. 10993-11	Standards Organization ANSI AAMI ISO	Standards Title Biological Evaluation of Medical Devices – Part 11. Tests for systemic toxicity	Version	Date 8/15/2006
5	Standards No. 11135-1	Standards Organization ANSI AAMI ISO	Standards Title Sterifization of Health Care Products Ethylene Oxide Part 1: Requirements for Development, Validation, and Routing Control of a Sterifization Process for Medical Devices	Version	Date 9/01/2007
6	Standards No. 11137-1	Standards Organization ANSI AAMI ISO	Standards Title Sterilization of Health Care Products – Radiation – Part 1; Requirements for development, validation and routine control of a sterilization process for medical devices	Version	Date 8/01/2006
7	Standards No.	Standards Organization ANSI AAMI ISO	Standards Title Sterilization of health care products - Radiation - Part 2: Establishing the sterilization dose	Version	2/15/2009

Please include any additional standards to be cited on a separate page.

Public reporting burden for this collection of information is estimated to average 0.5 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer 1350 Piccard Drive, Room 400 Rockville, MD 20850

An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number.

FORM FDA 3514 (12/10)

Page 5 of 5 Pages

EC		
_		

UTILIZATION OF STANDARDS

Note: Complete this section if your application or submission cites standards or includes a "Declaration of Conformity to a Recognized Standard" statement.

	Standards No.	Standards Organization	Standards Title	Version	Date
I	F136	ASTM	Standard Specification for Wrought Titanium-6Aluminum-4Vanadium ELI (Extra Low Interstitial) Alloy for Surgical Implant Applications		11/01/2008
	Standards No. F648	Standards Organization ASTM	Standards Title Standard Specification for Ultra-High-Molecular-Weight Polyethylene Powder and Fabricated Form for Surgical Implants	Version	Date 02/01/2007
2					02/01/2007
	Standards No.	Standards Organization	Standards Title Specifications for Wrought Cobalt-28 Chromium-6 Molybdenum	Version	Date
3		ASTM	Alloy for Surgical Implants		05/10/2000
4	Standards No.	Standards Organization	Standards Title	Version	Date
	Standards No.	Standards Organization	Standards Title	Version	Date
5					
6	Standards No.	Standards Organization	Standards Title	Version	Date
	Standards No.	Standards Organization	Standards Title	Version	Date
7					

Please include any additional standards to be cited on a separate page.

Public reporting hurden for this collection of information is estimated to average 0.5 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer 1350 Piccard Drive, Room 400 Rockville, MD 20850

An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number.

FORM FDA 3514 (12/10)

Page 5 of 5 Pages

Table of Contents

G	eneral and Administrative Information	1
	01. Medical Device User Fee Cover Sheet	1
	02. CDRH Premarket Review Submission Cover Sheet	2
Sι	ibmission	14
	03. Cover Letter	14
	04. Indications for Use Statement	15
	05. 510(k) Summary	16
	06. Truthful and Accurate Statement	18
	07. Class III Summary and Certification	19
	07a. Class III Summary	19
	Class III Summary and Certification	19
	08. Financial Certification/Disclosure	20
	08a. Financial Certification and Disclosure Statement	20
	Financial Certification or Disclosure Statement	20
	09. Declarations of Conformity and Summary Reports	21
	09a. Declaration of Conformity	21
	Declarations of Conformity and Summary Reports	21
	09b. Standards Data Report for 510(k)	22
	10. Executive Summary	56
	Device Name	56
	Reason for the Submission	56
	Device Description	56
	Overview	56
	Indications for Use	57
	Submission Overview	57
	Performance Testing	58
	i onormano i comg	20
	Stem Fatigue Testing	58
	-	
	Stem Fatigue Testing	58

Support Research	59
11. Device Description	61
11a. Device Description	61
Device Description	61
Overview	61
Humeral Component	63
Ulnar Component	64
Humeral Bearing-A	65
Ulnar Bearings-B	66
Axle Pin	68
Humeral Screws	69
System Compatibility	70
Size Interchangeability	70
Engineering Drawings	70
Catalog Numbers	71
Materials	72
Surface Characteristics	73
Mechanism of Action	74
Surgical Instrumentation Unique to the Device	75
Surgical Instrumentation Not Unique to the Device	83
Instrument Material Standards	84
Methods, Facilities and Controls	85
Method of Manufacturing	85
Packaging	86)
Humeral Stem	86
Ulnar Stem	86
Articulation Kit	87
Humeral Screw Kit	87
Pyrogenicity	87
Latex	87
11b. Engineering Drawings	88
12. Substantial Equivalence Discussion	114

3

Predicate Device(s)	114
Substantial Equivalence Comparison	114
Indications for Use	114
Technology	115
Performance Specifications	123
13. Proposed Labeling	130
13a. Labeling Summary	130
Labeling	130
13b. Draft Labels	131
13c. Draft Package Insert	136
13d. Draft Surgical Technique	143
Draft Surgical Technique	143
Pre-Operative Considerations	144
Short Version – Reference	146
Chapter I: Surgical Approach	147
Chapter 2: Humeral Preparation	149
Chapter 3: Ulnar Preparation	156
Chapter 4: Trial Reduction	160
Chapter 5: Bone Graft Preparation	160
Chapter 6: Ulnar Component Implantation	160
Chapter 7: Humeral Component Implantation	161
Chapter 8: Closure	166
Chapter 9: Postoperative Management	167
Chapter 10: Revision Technique	168
Chapter 11: Scrub Nurse Instructions	171
14. Sterilization and Shelf Life	176
Sterilization	176
Shelf Life	177
15. Biocompatibility	178
Biocompatibility	178
16. Software	179
Software	179

17. Electromagnetic Compatibility and Electrical Safety	180
Electromagnetic Compatibility and Electrical Safety	180
18. Performance Testing - Bench	181
18a. Bench Testing Summary	181
Performance Testing - Bench	181
Stem Fatigue Testing	182
Humeral Stem Fatigue Testing	182
Ulnar Stem, Mid-Stem Fatigue Testing	185
Ulnar Stem, Proximal Fatigue Testing	188
Wear Testing	191
Durability Testing	198
Modular Connection Fatigue Testing	205
Bench Testing in support of Vivacit-E Claim	207
18b. Test Reports	209
18b.1 ZRR_WA_2554_12	209
18b.2 ZRR_WA_2585_12	224
18b.3 ZRR_WA_2589_12	235
18b.4 ZRR_WA_2626_12	247
18b.5 ZRR_WA_2615_12	259
18b.6 ZRM_WA_0292_11	269
18b.7 ZRR_WA_2476_11	330
18b.8 TM1217.00	359
18b.9 TM1218.00	368
18b.10 ZRR_WA_2407_11	379
18b.11 ZRR_WA_2552_12	412
18b.12 ZRM_WA_0193_09	451
18b.13 ZRM_WA_0179_09	463
18b.14 ZRR_WA_1988_09	473
18b.15 ZRR_WA_2338_11	517
18b.16 TM 1222.00	532
18b.17 ZRM_WA_0300_11	537
18b.18 ZRR WA 2542 12	560

18b.19 ZRR_WI_2441_11	595
18b.20 ZRR_WA_2598_12	623
18b.21 ZRR_WA_2481_11	657
18b.22 Memo to DHF File Z07-014	668
18b.23 ZRR_WA_2252_10	670
18b.24 ZRR_WA_2648_12	685
18b.25 ZRR_WA_2646_12	705
18b.26 ZRR_WA_2580_12	723
18b.27 ZRR_WA_2409_11	744
18b.28 ZRR_WA_2399_11	77 2
18b.29 TM1267	794
18b.30 ZRR_WI_1222_12	800
19. Performance Testing - Animal	826
19a. Animal Testing Summary	826
Performance Testing - Animal	826
20. Performance Testing - Clinical	827
20a. Clinical Testing Summary	827
Performance Testing - Clinical	827
21. Letter of Authorization - MAF 1909	828
22. Letter of Authorization - MAF 652	829

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017



P.O. Box 708 Warsaw, IN 46581-0708 (574) 267-6131

13 December 2012

FDA CDRH DMC DEC 1 4 2012 Received

U.S. Food and Drug Administration Center for Devices and Radiological Health Document Mail Center – W066-G609 10903 New Hampshire Avenue Silver Spring, MD 20993-0002

Dear Sir or Madam:

Subject: Traditional 510(k) Premarket Notification – Zimmer Nexel Total Elbow – Prosthesis, Elbow, Constrained, Cemented – JDC – Class II – 21 CFR 888.3150

As required by Section 510(k) of the Federal Food, Drug and Cosmetic Act, as amended by the Medical Device Amendments of 1976, the Safe Medical Devices Act of 1990, FDA Modernization Act of 1997 (FDAMA) and in accordance with Title 21 of the Code of Federal Regulations (CFR) Part 807, subpart E, the above noted premarket notification is hereby submitted to the Food and Drug Administration (FDA). As required by 21 CFR 807.90(c), this document is submitted in duplicate. A CD-ROM containing an eCopy that is an exact duplicate of the paper copy submission is also included for your convenience.

These devices have not been previously submitted to FDA for identical or different indications, are not currently being reviewed for different indications by the same or different branch within ODE, and have not been previously cleared by FDA for different indications.

If you require any additional information or have any questions, please contact me by telephone at (574) 371-1642, e-mail at Joanna.Surma@Zimmer.com or fax at (574) 372-4605.

Sincerely,

form I durm 13 Dec 2012

Joanna L. Surma Associate Project Manager, Regulatory Affairs

Enclosure

Indications for Use

510(k) Number (if known):

Device Name:

Zimmer Nexel Total Elbow

Indications for Use:

Indications for use include:

- Elbow joint destruction which significantly compromises the activities of daily living
- Post-traumatic lesions or bone loss contributing to elbow instability
- Ankylosed joints, especially in cases of bilateral ankylosis from causes other than active sepsis
- Advanced rheumatoid, post-traumatic, or degenerative arthritis with incapacitating pain
- Instability or loss of motion when the degree of joint or soft tissue damage precludes reliable osteosynthesis
- Acute comminuted articular fracture of the elbow joint surfaces that precludes less radical procedures, including 13-C3 fractures of the distal humerus
- Revision arthroplasty

CAUTION: This device is intended for cemented use only.

Prescription Use <u>X</u> (Part 21 CFR 801 Subpart D) AND/OR

Over-The-Counter Use _____ (21 CFR 807 Subpart C)

(Please do not write below this line - Continue on another page if needed)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of 1

F



P.O. Box 708 Warsaw, IN 46581-0708 574 267-6131

Svh

Summary of Safety and Effectiveness

Sponsor:	Zimmer, Inc. P.O. Box 708 Warsaw, IN 46581-0708
Contact Person:	Joanna L. Surma Associate Project Manager, Regulatory Affairs Telephone: (574) 371-1642 Fax: (574) 372-4605
Date:	13 December 2012
Trade Name:	Zimmer Nexel Total Elbow
Common Name:	Total Elbow Prosthesis
Classification Name and Reference:	Prosthesis, Elbow, Constrained, Cemented 21 CFR § 888.3150
Product Code:	JDC
Predicate Device:	Coonrad/Morrey Total Elbow, manufactured by Zimmer, <u>K001989</u> , cleared 25 July 2000
	Coonrad/Morrey Total Elbow, manufactured by Zimmer, <u>K053189</u> , cleared 9 December 2005
Device Description:	The Zimmer Nexel total Elbow is a total elbow prosthesis designed for use with bone cement. It is available in multiple sizes and in right and left configurations.

Intended Use:

Indications for use include:

- Elbow joint destruction which significantly compromises the activities of daily living
- Post-traumatic lesions or bone loss contributing to elbow instability
- Ankylosed joints, especially in cases of bilateral ankylosis from causes other than active sepsis
- Advanced rheumatoid, post-traumatic, or degenerative arthritis with incapacitating pain
- Instability or loss of motion when the degree of joint or soft tissue damage precludes reliable osteosynthesis
- Acute comminuted articular fracture of the elbow joint surfaces that precludes less radical procedures, including 13-C3 fractures of the distal humerus
- Revision arthroplasty

Caution: This device is intended for cemented use only.

Page 2 of 2 13 December 2012

I.

Comparison to Predicate Device:	The Zimmer Nexel Total Elbow is substantially equivalent to the predicate device in terms of form and function. The Zimmer Nexel Total Elbow and the predicate device share similar intended uses and indications for use.
Performance Data:	 Non-Clinical Performance Testing Conducted: Stem Fatigue Testing Wear Testing Durability Testing Modular Connection Fatigue Testing
	Non-Clinical Performance Testing Conclusions:
	Non-clinical testing demonstrated that the New Zimmer Total Elbow meets performance requirements as defined by Design Control activities and is substantially equivalent to the predicate device in terms of safety and efficacy.

In this case, clinical data and conclusions were not needed to demonstrate substantial equivalence.



Traditional 510(k) Premarket Notification

TR

PREMARKET NOTIFICATION TRUTHFUL AND ACCURATE STATEMENT

[As Required by 21 CFR 807.87(k)]

Fertify that, in my capacity as <u>Regulatory Affairs Associate Project Manager</u> (title) of <u>Zimmer, Inc.</u> (company name), I believe to the best of my knowledge, that all data and information submitted in the premarket notification are truthful and accurate and that no material fact has been omitted.

12 Dec 2012 (Signature)

Joanna L. Surma (Typed Name)

13 Dec 2012 (Date)

(Premarket Notification [510(k)] Number)



Traditional 510(k) Premarket Notification

Class III Summary and Certification

This device is class II; therefore, this section does not apply.



Traditional 510(k) Premarket Notification

Financial Certification or Disclosure Statement

This submission does not contain any clinical data; therefore, this section does not apply.



1

Traditional 510(k) Premarket Notification

Declarations of Conformity and Summary Reports

This device conforms to the following recognized standards; a <u>Standards Data Report</u> (FDA Form 3654) is included below for each standard.

- ISO 14283 Implants for surgery Fundamental principles
- ISO 14630:2008 Non-active Surgical Implants General Requirements
- AAMI ST81:2004/(R)2010 Sterilization of medical devices Information to be provided by the manufacturer for the processing of resterilizable medical devices
- ANSI/AAMI/ISO 10993-1:2009 Biological Evaluation of Medical Devices Part 1: Evaluation and Testing within a risk management process
- ANSI/AAMI/ISO 10993-3:2003 Biological Evaluation of Medical Devices Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity
- ANSI/AAMI/ISO 10993-5:2009 Biological Evaluation of Medical Devices Part 5: Tests for in vitro cytotoxicity
- ANSI/AAMI/ISO 10993-6:2007 Biological Evaluation of Medical Devices Part 6: Tests for local effects after implantation
- ANSI/AAMI/ISO 10993-7:2008 Biological Evaluation of Medical Devices Part 7: Ethylene oxide sterilization residuals
- ANSI/AAMI/ISO 10993-10:2002/Amd.1:2006(E) Biological Evaluation of Medical Devices – Part 10: Tests for irritation and delayed-type hypersensitivity Amendment 1
- ANSI/AAMI/ISO 10993-11:2006 Biological Evaluation of Medical Devices Part 11: Tests for systemic toxicity
- ANSI/AAMI/ISO 11135-1:2007 Sterilization of Health Care Products Ethylene Oxide Part 1: Requirements for Development, Validation, and Routing Control of a Sterilization Process for Medical Devices
- ANSI/AAMI/ISO 11137-1:2006 Sterilization of Health Care Products Radiation – Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
- ANSI/AAMI/ISO 11137-2:2009 Sterilization of health care products Radiation – Part 2: Establishing the sterilization dose
- ASTM F136-08 Standard Specification for Wrought Titanium-6Aluminum-4Vanadium ELI (Extra Low Interstitial) Alloy for Surgical Implant Applications
- ASTM F648-07 Standard Specification for Ultra-High-Molecular-Weight Polyethylene Powder and Fabricated Form for Surgical Implants
- ASTM F1537-00 Specifications for Wrought Cobalt-28 Chromium-6 Molybdenum Alloy for Surgical Implants
- ASTM F1980-07 Standard Guidance for Accelerated Aging of Sterile Barrier Systems for Medical Devices

	Form Approved: OMB No. 0910-0120. Exp	iration Da	ate. 12/31/13
Food and Drug STANDARDS DATA	and Human Services Administration REPORT FOR 510(k)s n by applicant)		
This report and the Summary Report Table are to be comp ences a national or international standard. A separate repor			
TYPE OF 510(K) SUBMISSION			
STANDARD TITLE ' ISO 14283: 2004 Implants for Surgery – Fundamental Principles			
Please answer the following questions		Yes	No
Is this standard recognized by FDA ² ?			×
FDA Recognition number ³		ŧ	
Was a third party laboratory responsible for testing conform in the 510(k)?			X
Is a summary report ^a describing the extent of conformance 510(k)? If no, complete a summary report table.			×
Does the test data for this device demonstrate conformity to pertains to this device?		×	
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k)	······		X
Does this standard include more than one option or selection of selection of selected in the summary report table.	n of tests?		X
Were there any deviations or adaptations made in the use of If yes, were deviations in accordance with the FDA supplem			
Were deviations or adaptations made beyond what is specified of the specified of the section of the summar of the summar sections of adaptations in the summar sections of the section of			×
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			X
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation of Title of guidance:	dard? of this 510k?		
 The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication] Authonty [21 U S C 360d], http://www.fda.gov/MedicalDevices/ DeviceRegulationandGurdance/Standards/default htm http://www.accessdata.lda.gov/scnpts/cdrtvcfdocs/cfStandards/search.cfm 	address of the test laboratory or certification body invo assessment to this standard. The summary report incli- all standards utilized during the development of the de s The supplemental information sheet (SIS) is additional is necessary before FDA recognizes the standard. For www.accessdata Ida gov/scripts/cdrtvc/docs/cfStanda	udes infor wice Linformatic und at http	mation on on which o //
The summary report should include any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described, deviations from the standard, requirements not applicable to the device; and the name and	6 The online search for CDRH Guidance Documents cal http://www.fda.gov/Medica/Devices/DeviceRegulationa GuidanceDocuments/default htm	n be found	d al

FORM FDA 3654 (6/11)

-

ι

IN Belolug Network OF (2020) - FF

-		ANDARD CONFORMANCE RY REPORT TABLE		
STANDARD TITLE · ISO 14283: 2004 In	mplants for Surgery – Fundamental Prin	ciples		
<u></u>	CONFORMANCE V	VITH STANDARD SECTIONS'	<u> </u>	
SECTION NUMBER	SECTION TITLE		CONFORMANCE?	
All	Various		🕅 Yes 🔲 No	🗍 N/A
	OR OPTION SELECTED * to all indicated sections.		<u>,</u>	
DESCRIPTION				
JUSTIFICATION			···	
SECTION NUMBER	SECTION TITLE		CONFORMANCE?	
BEOTION THOMBEN				[_] N/A
TYPE OF DEVIATION	OR OPTION SELECTED *			
DESCRIPTION				
JUSTIFICATION				
SECTION NUMBER	SECTION	<u>. </u>	CONFORMANCE?	
TYPE OF DEVIATION	OR OPTION SELECTED *		Yes No	□ N/A
DESCRIPTION			<u></u>	
JUSTIFICATION				
explanation is need described and ade selected when follo report. More than * Types of deviations	ted under "justification." Some standards quately justified as appropriate for the su wing a standard is required under "type one page may be necessary. s can include an exclusion of a section in	ate whether conformance is met. If a section include options, so similar to deviations, t ibject device. Explanation of all deviations of deviation or option selected," "description in the standard, a deviation brought out by t to the device, or any adaptation of a section	the option chosen ne or description of option on" and "justification" the FDA supplementa	eds to be ons on the
	Paperwork R	teduction Act Statement		
time for revie completing at	wing instructions, searching existing dat	tion is estimated to average 1 hour per res ta sources, gathering and maintaining the on. Send comments regarding this burden gestions for reducing this burden to:	data needed, and	
Food Offic 1350	artment of Health and Human Services I and Drug Administration ee of Chief Information Officer I Piceard Drive, Room 400 kville, MD 20850	An agency may not conduct or spo required to respond to, a collection displays a currently valid OMB co	n of information unless	
ORM FDA 3654 (6/1	1)	Page 2		

1

	Form Approved: OMB No. 0910-0120; Ex	piration Da	te: 12/31/13
Food and Drug STANDARDS DATA F	and Human Services Administration REPORT FOR 510(k)s n by applicant)		
This report and the Summary Report Table are to be comp ences a national or international standard. A separate report			
TYPE OF 510(K) SUBMISSION	Abbreviated		
STANDARD TITLE ' ISO 14630: 2008 Non-active Surgical Implants - General Require	ements		
Please answer the following questions		Yes	No
Is this standard recognized by FDA 2?		×	
FDA Recognition number ³		# <u>11-208</u>	
Was a third party laboratory responsible for testing conformining the 510(k)?			×
Is a summary report ⁴ describing the extent of conformance 510(k)? If no, complete a summary report table.			X
Does the test data for this device demonstrate conformity to pertains to this device?		R	
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).			×
Does this standard include more than one option or selection If yes, report options selected in the summary report table.	n of tests?		X
Were there any deviations or adaptations made in the use o If yes, were deviations in accordance with the FDA supplem			
Were deviations or adaptations made beyond what is specif If yes, report these deviations or adaptations in the summary			X
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			X
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation o Title of guidance:	of this 510k?		
 ¹ The formatting convention for the title is [SDO] [numeric identifier] [title of standard] [date of publication] ² Authority [21 U.S.C. 360d], http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm ³ http://www.accessdata.fda.gov/scnpts/cdrh/ctdocs/cfStandards/search.cfm ⁴ The summary report should include, any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described, deviations from the standard, requirements not applicable to the device, and the name and 	address of the test laboratory or certification body inv assessment to this standard. The summary report ind all standards utilized during the development of the d s The supplemental information sheet (SIS) is addition is necessary before FDA recognizes the standard. Fo www.accessdata.fda.gov/scnpts/cdrh/cfdocs/cfStand e The online search for CDRH Guidance Documents ca http://www.fda.gov/Medica/Devices/DeviceRegulation GuidanceDocuments/default htm	cludes inform levice, al informatio bund at http: ards/search an be found	mation on on which :// i.cfm

FORM FDA 3654 (6/11)

ц.

x.

The tradition ${\bf S}_{\rm eff}$, we show that the ${\bf M}_{\rm eff}={\bf F}$ (

		NDARD CONFORMANCE Y REPORT TABLE
STANDARD TITLE ISO 14630: 2008 N	Ion-active Surgical Implants – General R	equirements
	CONFORMANCE W	TH STANDARD SECTIONS
SECTION NUMBER	SECTION TITLE	CONFORMANCE?
All	Various	🔀 Yes 🗌 No 🗍 N/A
	OR OPTION SELECTED * to all indicated sections.	
DESCRIPTION		
JUSTIFICATION		
SECTION NUMBER	SECTION TITLE	CONFORMANCE?
		Yes No N/A
TYPE OF DEVIATION	OR OPTION SELECTED *	
DESCRIPTION		
JUSTIFICATION		
SECTION NUMBER	SECTION TITLE	CONFORMANCE?
4		□ Yes □ N₀ □ N/A
TYPE OF DEVIATION	OR OPTION SELECTED *	
DESCRIPTION		
JUSTIFICATION	···-	······
explanation is need described and adec selected when follo report. More than c * Types of deviations	led under "justification " Some standards quately justified as appropriate for the sub wing a standard is required under "type o one page may be necessary. s can include an exclusion of a section in	te whether conformance is met. If a section is not applicable (N/A) an include options, so similar to deviations, the option chosen needs to be opect device. Explanation of all deviations or description of options if deviation or option selected," "description" and "justification" on the the standard, a deviation brought out by the FDA supplemental the device, or any adaptation of a section.
	Paperwork R	eduction Act Statement
time for revie completing an	wing instructions, searching existing data	on is estimated to average 1 hour per response, including the a sources, gathering and maintaining the data needed, and n. Send comments regarding this burden estimate or any other restions for reducing this burden to:
Depc	artment of Health and Human Services	
	and Drug Administration	
	e of Chief Information Officer Piccard Drive, Room 400	An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it
	ville, MD 20850	displays a currently valid OMB control number.

i

-

	Form Approved: OMB No. 0910-0120: Exp	piration Da	te: 12/31/13
Food and Drug STANDARDS DATA	and Human Services g Administration REPORT FOR 510(k)s n by applicant)		
This report and the Summary Report Table are to be comp ences a national or international standard. A separate repor			
TYPE OF 510(K) SUBMISSION			,
STANDARD TITLE ' ASTM F1980-07 Standard Guidance for Accelerated Aging of Ster	rile Barrier Systems for Medical Devices		
Please answer the following questions		Yes	No
Is this standard recognized by FDA ??			
FDA Recognition number 3	· · · · · · · · · · · · · · · · · · ·	¥14-229	
Was a third party laboratory responsible for testing conforming in the 510(k)?			X
Is a summary report ⁴ describing the extent of conformance 510(k)? If no, complete a summary report table.			
Does the test data for this device demonstrate conformity to pertains to this device?		X	
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).			$\mathbf{\Sigma}$
Does this standard include more than one option or selectio If yes, report options selected in the summary report table.	on of tests?		X
Were there any deviations or adaptations made in the use of If yes, were deviations in accordance with the FDA supplement			
Were deviations or adaptations made beyond what is specif If yes, report these deviations or adaptations in the summar			X
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation of Title of guidance:	of this 510k?		
 The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication] Authonty [21 U.S.C. 360d], http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm The summary report should include, any adaptations used to adapt to the 	address of the test laboratory or certification body invo assessment to this standard. The summary report incl all standards utilized during the development of the de s The supplemental information sheet (SIS) is additional is necessary before FDA recognizes the standard. Fo www.accessdata fda gov/scnpts/cdrh/cfdocs/cfStanda a The online search for CDRH Guidance Documents ca	ludes inform avice. Il informatio und at http: irds/search	mation on on which ;// i cfm
device under review (for example, alternative test methods): choices made when options or a selection of methods are described; deviations from the standard, requirements not applicable to the device, and the name and	http://www.fda.gov/MedicalDevices/DeviceRegulation GuidanceDocuments/default.htm	andGuidan	ce/

FORM FDA 3654 (6/11)

ĵ,

USC Bitlistang Services (2017) 1977 (2017)

1

		NDARD CONFORMANCE Y REPORT TABLE	
STANDARD TITLE ASTM F1980-07 Sta	andard Guidance for Accelerated Aging o	of Sterile Barrier Systems for Medical Dev	ices
	CONFORMANCE W	ITH STANDARD SECTIONS*	
SECTION NUMBER	SECTION TITLE		CONFORMANCE?
All	Various		🕅 Yes 🗌 No 📋 N/A
	OR OPTION SELECTED * to all indicated sections.		
DESCRIPTION		·····	
JUSTIFICATION			
SECTION NUMBER	SECTION TITLE		CONFORMANCE?
			Yes No N/A
TYPE OF DEVIATION	OR OPTION SELECTED *		
DESCRIPTION			
JUSTIFICATION			
SECTION NUMBER	SECTION TITLE	· · ·	CONFORMANCE?
			Yes No N/A
TYPE OF DEVIATION	OR OPTION SELECTED *		
DESCRIPTION			
JUSTIFICATION			
explanation is need described and aded selected when follo report. More than d * Types of deviations	led under "justification." Some standards i quately justified as appropriate for the sub wing a standard is required under "type o one page may be necessary. s can include an exclusion of a section in	te whether conformance is met. If a section include options, so similar to deviations, th oject device. Explanation of all deviations o if deviation or option selected." "description the standard, a deviation brought out by th	e option chosen needs to be r description of options " and "justification" on the
information sheet (3	SIS), a deviation to adapt the standard to	the device. or any adaptation of a section.	
	•	duction Act Statement	
time for revie completing ar	wing instructions, searching existing data	on is estimated to average 1 hour per resp a sources, gathering and maintaining the da n. Send comments regarding this burden es- sections for reducing this burden to:	ita needed, and
Food Offic 1350	artment of Health and Human Services I and Drug Administration are of Chief Information Officer Piccard Drive, Room 400 wille, MD 20850	An agency may not conduct or spon required to respond to, a collection displays a currently valid OMB con	of information unless it

.

Form Approved: OMB No. 0910-0120; Ex	xpiration Date, 12/31	1/13
--------------------------------------	-----------------------	------

Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)					
This report and the Summary Report Table are to be com ences a national or international standard. A separate repo	pleted by the applicant when submitting a rt is required for each standard referenced	510(k) t in the 51	hat refer- 10(k).		
TYPE OF 510(K) SUBMISSION	n				
Traditional Special	Abbreviated				
STANDARD TITLE ' ANSI/AAMI ST81:2004/(R)2010 Sterilization of medical devices		irer for t	he proc		
Please answer the following questions		Yes	No		
Is this standard recognized by FDA 2?		X			
FDA Recognition number ³		¥14-295			
Was a third party laboratory responsible for testing conform in the 510(k)?			Ø		
Is a summary report ⁴ describing the extent of conformance 510(k)? If no, complete a summary report table.					
Does the test data for this device demonstrate conformity t pertains to this device?			×		
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).			X		
Does this standard include more than one option or selection of selection of selected in the summary report table.					
Were there any deviations or adaptations made in the use If yes, were deviations in accordance with the FDA supplement					
Were deviations or adaptations made beyond what is spec If yes, report these deviations or adaptations in the summa			X		
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			X		
Is there an FDA guidance ^s that is associated with this stan If yes, was the guidance document followed in preparation Title of guidance:	of this 510k?				
 ¹ The formatting convention for the title is [SDO] (numeric identifier) [title of standard) [date of publication] ² Authority [21 U.S.C. 360d], http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm 	address of the test laboratory or certification body inve assessment to this standard. The summary report inc all standards utilized during the development of the de s The supplemental information sheet (SIS) is additional is necessary before FDA recognizes the standard. Fo	ludes infor avice Il informati	malion on on which		
³ http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm ⁴ The summary report should include, any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device, and the name and	Is necessary bacine FDX recognizes the standard. Fo www accessidata. Ida.gov/scripts/cdth/cfdocs/cfStandard. The online search for CDRH Guidance Documents ca http://www.fda.gov/MedicalDevices/DeviceRegulation GuidanceDocuments/default.htm	irds/searcl	h.cfm diat		
FORM FDA 3654 (6/11)	ge 1	iblinhing Service	nywytatatau - E		

ŧ

EXTENT OF STANDARD CONFORMANCE
SUMMARY REPORT TABLE

					<u> </u>
STANDARD TITLE ANSI/AAMI ST81:20	04/(R)2010 Sterilization of medical dev	vices — Information to be provided by the	manufactu	rer for the	е ргос
	CONFORMANCE W	ITH STANDARD SECTIONS*			
SECTION NUMBER	SECTION TITLE	·····	CONFORM	ANCE?	
All Sections	Various		🔀 Yes	No No	[] N/A
	R OPTION SELECTED * all indicated sections.				
DESCRIPTION					
JUSTIFICATION	·····				
SECTION NUMBER	SECTION TITLE		CONFORM	ANCE?	
			🗌 Yes	No No	🗌 N/A
TYPE OF DEVIATION O	R OPTION SELECTED +		··		
DESCRIPTION					
JUSTIFICATION					
SECTION NUMBER	SECTION TITLE		CONFORM	ANCE?	
			🗌 Yes	No No	🗌 N/A
TYPE OF DEVIATION O	R OPTION SELECTED +		. <u>.</u>		
DESCRIPTION					
JUSTIFICATION					
explanation is neede described and adequ selected when follow report. More than on * Types of deviations of	d under "justification." Some standards i lately justified as appropriate for the sub ring a standard is required under "type o le page may be necessary. can include an exclusion of a section in	e whether conformance is met. If a sectior include options, so similar to deviations, th ject device. Explanation of all deviations of f deviation or option selected," "description the standard, a deviation brought out by th the device, or any adaptation of a section.	e option ch or descriptio n" and "justi ae FDA supp	osen nee n of optic fication" o	eds to be ons on the
	Paperwork Re	duction Act Statement			
time for review completing and	g burden for this collection of informati ing instructions, searching existing data	on is estimated to average 1 hour per resp a sources, gathering and maintaining the d n. Send comments regarding this burden e	ata needed,	and	
Food a Office 1350 I	tment of Health and Human Services and Drug Administration of Chief Information Officer Piccard Drive, Room 400 rille, MD 20850	An agency may not conduct or spor required to respond to, a collection displays a currently valid OMB con	of informati	on unless	
EORM EDA 3654 (6/11)	······································	Page 2			

ï

Form Approved: OMB No. 0910-0120; Expiration Date: 12/31/13

30

Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)					
This report and the Summary Report Table are to be comp ences a national or international standard. A separate report					
TYPE OF 510(K) SUBMISSION					
🗙 Traditionał 🔄 Special	Abbrevialed				
STANDARD TITLE ' ANSI/AAMI/ISO 10993-1:2009 Biological Evaluation of Medical	•	a Risk i	Manage		
Please answer the following questions		Yes	No		
Is this standard recognized by FDA 2?		\boxtimes			
FDA Recognition number ³		<u>‡2-156</u>			
Was a third party laboratory responsible for testing conform in the 510(k)?					
Is a summary report ⁴ describing the extent of conformance 510(k)? If no, complete a summary report table.			×		
Does the test data for this device demonstrate conformity to pertains to this device?		X			
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).			×.		
Does this standard include more than one option or selection If yes, report options selected in the summary report table.	on of tests?	×			
Were there any deviations or adaptations made in the use of If yes, were deviations in accordance with the FDA supplement					
Were deviations or adaptations made beyond what is speci If yes, report these deviations or adaptations in the summar			X		
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			×		
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation	of this 510k?				
Title of guidance:					
 The formatting convention for the title is: [SDO] [numenc identifier] [title of standard] [date of publication] 	address of the test laboratory or certification body invo assessment to this standard. The summary report incl all standards utilized during the development of the de	udes infor			
² Authority [21 U S.C. 360d], http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm 3 http://www.saccodolo.fdp.com/coccod/default.htm	The supplemental information sheet (SIS) is additional is necessary before FDA recognizes the standard. For	und at http	.H		
 ³ http://www.accessdata.fda.gov/scnpts/cdrh/cfdocs/cfStandards/search.cfm ⁴ The summary report should include, any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and 	www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStanda e The online search for CDRH Guidance Documents ca http://www.fda.gov/MedicalDevices/DeviceRegulationa GuidanceDocuments/default.htm	n be found	l at		
FORM FDA 3654 (6/11) Pag	је 1 вск	bhahag Services	and the El		

t

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE				
STANDARD TITLE ANSI/AAMI/ISO 109	93-1:2009 Biological Evaluation of M	fedical Devices - Part 1: Evaluation and Te	sting within a Risk Manage	
	CONFORMANCE	WITH STANDARD SECTIONS*		
SECTION NUMBER All Sections	SECTION TITLE Various		CONFORMANCE?	
A biological evaluation	R OPTION SELECTED * on was conducted per ISO 10993-1, wi nsitization, 5.2.5 Intracutaneous React	hich led to the selected ISO 10993 testing o tivity, 5.2.6 Systemic Toxicity		
JUSTIFICATION				
SECTION NUMBER	SECTION TITLE			
TYPE OF DEVIATION O	R OPTION SELECTED *			
DESCRIPTION		· · · · · · · · · · · · · · · · · · ·		
JUSTIFICATION				
SECTION NUMBER	SECTION TITLE		CONFORMANCE?	
TYPE OF DEVIATION O	R OPTION SELECTED *		Yes No N/A	
DESCRIPTION	,			
JUSTIFICATION				
explanation is needed described and adequiselected when follow report. More than or * Types of deviations	ed under "justification." Some standard uately justified as appropriate for the s ving a standard is required under "type ne page may be necessary. can include an exclusion of a section i	ate whether conformance is met. If a section is include options, so similar to deviations, the ubject device. Explanation of all deviations of of deviation or option selected," "description in the standard, a deviation brought out by the to the device, or any adaptation of a section	ne option chosen needs to be or description of options n" and "justification" on the ne FDA supplemental	
	Paperwork I	Reduction Act Statement		
Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to:				
Food Office 1350	rtment of Health and Human Services and Drug Administration e of Chief Information Officer Piccard Drive, Room 400 ville, MD 20850	An agency may not conduct or spor required to respond to, a collection displays a currently valid OMB cor	of information unless it	

Carrier Assessments OHO Mar	AAAA AAAA Custoniise Deter	40104140
FORM ADDROVED CIMIN NO.	0910-0120; Expiration Date:	120.1101
i onn approved. Onio no.	borto orgo, Expiration Date.	

Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)					
This report and the Summary Report Table are to be comp ences a national or international standard. A separate report					
TYPE OF 510(K) SUBMISSION					
Traditional Special	Abbreviated				
STANDARD TITLE ' ANSI/AAMI/ISO 10993-3:2003/(R)2009 Biological Evaluation of	Medical Devices - Part 3: Tests for Genotoxici	ty, Carci	inogenic		
Please answer the following questions		Yes	No		
Is this standard recognized by FDA 2?		×			
FDA Recognition number ³	······································	<u><u></u>¥2-117</u>			
Was a third party laboratory responsible for testing conforming the 510(k)?	1. Andres allows and Analysis — Statistic control and the second control of the second control and the second control of the seco	X			
Is a summary report ⁴ describing the extent of conformance 510(k)? If no, complete a summary report table.			×		
Does the test data for this device demonstrate conformity to pertains to this device?		X			
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).	······	×			
Does this standard include more than one option or selection of selection of the summary report table.	n of tests?	×			
Were there any deviations or adaptations made in the use of	f the standard?				
If yes, were deviations in accordance with the FDA supplem	iental information sheet (SIS) 5?		\Box		
Were deviations or adaptations made beyond what is specifing of the section of the section of the summar section of the sectio			\boxtimes		
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.					
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation of Title of guidance:	of this 510k?				
¹ The formatting convention for the title is. [SDO] [numeric identifier] [title of standard] [date of publication]	address of the test laboratory or certification body invo assessment to this standard. The summary report incl	udes infor			
² Authority [21 U.S.C. 360d], http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm	all standards utilized during the development of the de • The supplemental information sheet (SIS) is additional	l informati			
³ http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm	is necessary before FDA recognizes the standard. For www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStanda				
The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard, requirements not applicable to the device; and the name and	6 The online search for CDRH Guidance Documents ca http://www.fda.gov/Medica/Devices/DeviceRegulation. GuidanceDocuments/default.htm				
FORM FDA 3654 (6/11) Pag	e1 rscit	ilahay serve	շարմությու ը		
Questions? Contact FDA/CDRH/OCE/DID at CDF	RH-FOISTATUS@fda.hhs.gov or 301-796	-8118	32		

۰.

ŕ

EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE					
STANDARD TITLE ANSI/AAMI/ISO 109	93-3:2003/(R)2009 Biological Evaluation	of Medical Devices - Part 3: Tests for (Genotoxicit	y, Carcir	ogenic
	CONFORMANCE WITH	STANDARD SECTIONS			
SECTION NUMBER All Sections	SECTION TITLE Various			ANCE?	
TYPE OF DEVIATION OF Test options selected in Study DESCRIPTION	ROPTION SELECTED * nclude: Bacterial Reverse Mutation, Mous	e Lymphoma Assay, and Mouse Perip!		<u> </u>	
JUSTIFICATION					
SECTION NUMBER	SECTION TITLE			IANCE?	□ N/A
TYPE OF DEVIATION OF	R OPTION SELECTED +			<u> </u>	
DESCRIPTION					
JUSTIFICATION					
SECTION NUMBER	SECTION TITLE				
TYPE OF DEVIATION O	R OPTION SELECTED +		Yes	No	N/A
DESCRIPTION					
JUSTIFICATION					
 * For completeness list all sections of the standard and indicate whether conformance is met. If a section is not applicable (N/A) an explanation is needed under "justification." Some standards include options, so similar to deviations, the option chosen needs to be described and adequately justified as appropriate for the subject device. Explanation of all deviations or description of options selected when following a standard is required under "type of deviation or option selected," "description" and "justification" on the report. More than one page may be necessary. * Types of deviations can include an exclusion of a section in the standard, a deviation brought out by the FDA supplemental information sheet (SIS), a deviation to adapt the standard to the device, or any adaptation of a section. 					
	Paperwork Redu	ction Act Statement			
time for review completing and	Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to:				
Food a Office 1350 F	tment of Health and Human Services and Drug Administration of Chief Information Officer Piccard Drive, Room 400 ille, MD 20850	An agency may not conduct or spon required to respond to, a collection displays a currently valid OMB con	of informatio		

ŝ,

Form Approved: OMB No. 0910-0120; Expiration Date: 12/31/13

Food and Drug STANDARDS DATA	n and Human Services g Administration REPORT FOR 510(k)s n by applicant)		
This report and the Summary Report Table are to be comp ences a national or international standard. A separate report			
TYPE OF 510(K) SUBMISSION			_
Traditional Special			_
STANDARD TITLE 1 ANSI/AAMI/ISO 10993-5:2009 Biological Evaluation of Medical	Devices - Part 5: Tests for In Vitro Cytotoxicit	у	
Please answer the following questions		Yes	No
Is this standard recognized by FDA 2?		×	
FDA Recognition number ³		¥2-153	
Was a third party laboratory responsible for testing conform in the 510(k)?			
Is a summary report ⁴ describing the extent of conformance 510(k)? If no, complete a summary report table.			×
Does the test data for this device demonstrate conformity to pertains to this device?	433 (1973) 10713 10713 103 (193, 193, 193, 193, 193, 193, 193, 193,	×	
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).		\boxtimes	
Does this standard include more than one option or selection If yes, report options selected in the summary report table.	on of tests?	X	
Were there any deviations or adaptations made in the use of If yes, were deviations in accordance with the FDA supplement			
Were deviations or adaptations made beyond what is specified of the summar lift yes, report these deviations or adaptations in the summar			
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			×
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation of Title of guidance:	of this 510k?		
 ¹ The formatting convention for the title is: (SDO) [numeric identifier] [title of standard] [date of publication] ² Authority [21 U S.C. 360d], http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default htm ³ http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm ⁴ The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods), choices made when options or a selection of methods are described; daviations from the standard, requirements not applicable to the device, and the name and 	address of the test laboratory or certification body invi assessment to this standard. The summary report inc all standards utilized during the development of the de 5 The supplemental information sheet (SIS) is additional is necessary before FDA recognizes the standard. Fo www.accessdata.fda.gov/scripts/cdm/cfdocs/cfStanda 6 The online search for CDRH Guidance Documents ca http://www.fda.gov/MedicalDevices/DeviceRegulation GuidanceDocuments/default.htm	ludes infor avice. Il information und at http ards/search in be found	mation on on which y II n.cfm d at

L

Page 1

	EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE					
STANDARD TITLE ANSI/AAMI/ISO 109	93-5:2009 Biological Evaluation of Med	lical Devices - Part 5: Tests for In Vitro (
	CONFORMANCE WI	TH STANDARD SECTIONS"				
SECTION NUMBER All Sections	SECTION TITLE Various					
	I R OPTION SELECTED • nclude a cytotoxicity study using the elu	tion method				
DESCRIPTION						
JUSTIFICATION						
SECTION NUMBER	SECTION TITLE					
TYPE OF DEVIATION O	R OPTION SELECTED +	······				
DESCRIPTION						
JUSTIFICATION						
SECTION NUMBER	SECTION TITLE	<u> </u>	CONFORMANCE?			
TYPE OF DEVIATION O	R OPTION SELECTED +		_ Yes _ No _ N/A			
DESCRIPTION						
JUSTIFICATION						
explanation is neede described and adequ selected when follow report. More than or * Types of deviations of	d under "justification." Some standards in lately justified as appropriate for the subj ing a standard is required under "type of le page may be necessary. can include an exclusion of a section in th	whether conformance is met. If a section include options, so similar to deviations, th act device. Explanation of all deviations o deviation or option selected," "description ne standard, a deviation brought out by th he device, or any adaptation of a section.	e option chosen needs to be or description of options " and "justification" on the e FDA supplemental			
Paperwork Reduction Act Statement Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to:						
Food a Office 1350 I	tment of Health and Human Services and Drug Administration of Chief Information Officer Piccard Drive, Room 400 ville, MD 20850	An agency may not conduct or spon required to respond to, a collection displays a currently valid OMB con	of information unless it			

ī

Form Approved: OMB No. 0910-0120; I	Expiration Da	ate: <u>12/3</u> 1/13
Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)		
This report and the Summary Report Table are to be completed by the applicant when submitting ences a national or international standard. A separate report is required for each standard reference		
TYPE OF 510(K) SUBMISSION		
STANDARD TITLE 1 ANSI/AAMI/ISO 10993-6:2007 Biological Evaluation of Medical Devices - Part 6: Tests for Local Effects after	er Implanta	tion
Please answer the following questions	Yes	No
Is this standard recognized by FDA 2?		
FDA Recognition number ³	#2-120	
Was a third party laboratory responsible for testing conformity of the device to this standard identifie in the 510(k)?		
Is a summary report ⁴ describing the extent of conformance of the standard used included in the 510(k)? If no, complete a summary report table.		×
Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device?		
Does this standard include acceptance criteria?	×	
Does this standard include more than one option or selection of tests? If yes, report options selected in the summary report table.		
Were there any deviations or adaptations made in the use of the standard? If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) ⁵ ?	· · · · · ·	
Were deviations or adaptations made beyond what is specified in the FDA SIS? If yes, report these deviations or adaptations in the summary report table.	· 🔲	×
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.		×
Is there an FDA guidance ⁶ that is associated with this standard? If yes, was the guidance document followed in preparation of this 510k? Title of guidance:		
 The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication] Authority [21 U S.C 360d], http://www fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm The summary report should include, any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard, requirements not applicable to the device, and the name and address of the test laboratory or certification body assessment to this standard. The summary report all standards utilized during the development of the standard www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm The summary report should include, any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard, requirements not applicable to the device, and the name and 	includes infor e device onal informati Found at http ndards/searc s can be found	mation on ion which o.// h.cfm d at

FORM FDA 3654 (6/11)

Page 1

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118 36

EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE				
STANDARD TITLE ANSI/AAMI/ISO 109	93-6:2007 Biological Evaluation of M	ledical Devices - Part 6: Tests for Local Eff	ects after Implantation	
	CONFORMANCE	MITH STANDARD SECTIONS*		
SECTION NUMBER	SECTION TITLE		CONFORMANCE?	
	PR OPTION SELECTED *		Yes No N	NA
	include muscle implantation studies in	rabbits at 2 and 12 weeks		
DESCRIPTION	<u> </u>			
JUSTIFICATION				
SECTION NUMBER	SECTION TITLE		CONFORMANCE?	
	_		Yes No N	N/A
TYPE OF DEVIATION O	R OPTION SELECTED *			
DESCRIPTION				
JUSTIFICATION				
SECTION NUMBER	SECTION TITLE		CONFORMANCE?	
	R OPTION SELECTED +			1/A
TYPE OF DEVIATION C	R OPTION SELECTED			i
DESCRIPTION	, <u></u>			
JUSTIFICATION	<u> </u>			
 For completeness list all sections of the standard and indicate whether conformance is met. If a section is not applicable (N/A) an explanation is needed under "justification." Some standards include options, so similar to deviations, the option chosen needs to be described and adequately justified as appropriate for the subject device. Explanation of all deviations or description of options selected when following a standard is required under "type of deviation or option selected," "description" and "justification" on the report. More than one page may be necessary. Types of deviations can include an exclusion of a section in the standard, a deviation brought out by the FDA supplemental information sheet (SIS), a deviation to adapt the standard to the device, or any adaptation of a section 				
	Paperwork P	Reduction Act Statement		_
Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to:				
Food Office 1350	rtment of Health and Human Services and Drug Administration e of Chief Information Officer Piccard Drive, Room 400 ville, MD 20850	An agency may not conduct or spon required to respond to, a collection displays a currently valid OMB con	of information unless it	
				_

Form Approved: OMB No. 0910-0120; Expiration Date: 12/31/13					
Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)					
This report and the Summary Report Table are to be completed by the applicant when submitting a ences a national or international standard. A separate report is required for each standard referenced is					
TYPE OF 510(K) SUBMISSION					
STANDARD TITLE 1 ANSI/AAMI/ISO 10993-7:2008 Biological Evaluation of Medical Devices - Part 7: Ethylene Oxide Sterilization	Residuals	;			
Please answer the following questions	Yes	No			
Is this standard recognized by FDA ² ?	X				
FDA Recognition number ³	⊭ 14-278				
Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)?	×				
Is a summary report ⁴ describing the extent of conformance of the standard used included in the 510(k)? If no, complete a summary report table.					
Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device?	×				
Does this standard include acceptance criteria?	×				
Does this standard include more than one option or selection of tests? If yes, report options selected in the summary report table.	X				
Were there any deviations or adaptations made in the use of the standard? If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) ⁵ ?					
Were deviations or adaptations made beyond what is specified in the FDA SIS? If yes, report these deviations or adaptations in the summary report table.		X			
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.		×			
Is there an FDA guidance ⁶ that is associated with this standard? If yes, was the guidance document followed in preparation of this 510k? Title of guidance:					
 ¹ The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication] ² Authority [21 U.S C. 360d]. http://www fda gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm ³ http://www accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search cfm ⁴ The summary report should include. any adaptations used to adapt to the device under review (for example, alternative test methods); choices made ⁵ The online search for CDRH Guidance Documents can be online search for CDRH Guidance Documents	ludes inform evice Il informatio und at http: ards/search. In be found	nation on m which // .cfm at			
when options or a selection of methods are described; deviations from the GurdanceDocuments/default htm	andGuidani	ce/			

standard; requirements not applicable to the device, and the name and

1

Page 1

INCINHABAN SERVICE FOR LARGE OF EF

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118 38

GuidanceDocuments/default.htm

EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE						
STANDARD TITLE ANSI/AAMI/ISO 109	STANDARD TITLE ANSI/AAMI/ISO 10993-7:2008 Biological Evaluation of Medical Devices - Part 7: Ethylene Oxide Sterilization Residuals					
	CONFORMANCE WITH	STANDARD SECTIONS*				
SECTION NUMBER	SECTION TITLE		CONFORMANCE?			
4.2 & 4.3	4.2) Categorization of Medical Devices 4	4.3) Allowable Limits	X Yes No N/A			
TYPE OF DEVIATION OF The Vitamin E HXPE contact).	R OPTION SELECTED * Bearings and the Zimaloy Axle-Pin are ca	legorized as permanent contact devices	(greater than 30 days of			
DESCRIPTION The EO and ECH limi	ts for permanent contact devices were utili	zed				
JUSTIFICATION						
SECTION NUMBER	SECTION TITLE		CONFORMANCE?			
4.4	Determination of EO and ECH Residuals		Yes 🗌 No 🗍 N/A			
TYPE OF DEVIATION O Exhaustive Extraction	R OPTION SELECTED * Method is utilized instead of simulated us	e method				
	e limits are met & residues shown to be wi	thin limits for exhaustive extraction the	ere is no need to test with			
simulated use						
3031110/11011						
SECTION NUMBER	SECTION TITLE		CONFORMANCE?			
All others	Various		Yes No N/A			
TYPE OF DEVIATION O N/A	R OPTION SELECTED *					
DESCRIPTION None, Conformance to	all indicated sections.					
JUSTIFICATION						
 * For completeness list all sections of the standard and indicate whether conformance is met. If a section is not applicable (N/A) an explanation is needed under "justification." Some standards include options, so similar to deviations, the option chosen needs to be described and adequately justified as appropriate for the subject device. Explanation of all deviations or description of options selected when following a standard is required under "type of deviation or option selected," "description" and "justification" on the report. More than one page may be necessary * Types of deviations can include an exclusion of a section in the standard, a deviation brought out by the FDA supplemental information sheet (SIS), a deviation to adapt the standard to the device, or any adaptation of a section. 						
	Paperwork Redu	iction Act Statement				
Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to:						
Food a Office 1350 (aspect of this collection of information, including suggestions for reducing this burden to: Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer 1350 Piccard Drive, Room 400 Rockville, MD 20850 An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number.					

Is this standard recognized by FDA ² ? Image: Complete a summary report table. FDA Recognition number ³ #2-152 Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? Image: Complete a summary report 4 describing the extent of conformance of the standard used included in the 510(k)? If no, complete a summary report table. Image: Complete a standard as it	k).
ences a national or international standard. A separate report is required for each standard referenced in the 510(TYPE OF 510(K) SUBMISSION Image: Traditional Image: Traditional STANDARD TITLE ' ISO 10993-10:2002/Amd. 1:2006(E) Biological Evaluation of Medical Devices - Part 10: Tests for Irritation and Delayed-Ty Please answer the following questions Yes Is this standard recognized by FDA ² ? Image: Test for Irritation and Delayed (Image: Test for Irritation (Image: Test for Irritati	k).
Image: Special in the standard used included in the standard used included in the story report table. Special in the standard used included in the standard as it	pe
ISO 10993-10:2002/Amd. 1:2006(E) Biological Evaluation of Medical Devices - Part 10: Tests for Irritation and Delayed-Ty Please answer the following questions Yes Is this standard recognized by FDA ² ? Image: Standard recognized by FDA ² ? FDA Recognition number ³ #2-152 Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? Image: Standard used included in the 510(k)? Is a summary report ⁴ describing the extent of conformance of the standard used included in the 510(k)? Image: Standard used included in the 510(k)? If no, complete a summary report table. Image: Standard used included in the 510(k)? Does the test data for this device demonstrate conformity to the requirements of this standard as it	pe
Is this standard recognized by FDA ² ? Image: Standard recognized by FDA ² ? FDA Recognition number ³ #2-152 Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? Image: Standard Identified Image: Standard	
FDA Recognition number 3 #2-152 Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? Is a summary report 4 describing the extent of conformance of the standard used included in the 510(k)? Is a summary report 4 describing the extent of conformance of the standard used included in the 510(k)? Image: Conformance of the standard used included in the 510(k)? If no, complete a summary report table. Image: Conformation of the standard as it	No
Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? Image: Conformance of the standard used included in the 510(k)? Is a summary report ⁴ describing the extent of conformance of the standard used included in the 510(k)? Image: Conformance of the standard used included in the 510(k)? If no, complete a summary report table. Image: Conformance of the standard used included in the 510(k)? Does the test data for this device demonstrate conformity to the requirements of this standard as it	
in the 510(k)? Is a summary report ⁴ describing the extent of conformance of the standard used included in the 510(k)? If no, complete a summary report table.	
510(k)?	
	X
Does this standard include acceptance criteria?	
Does this standard include more than one option or selection of tests?	
Were there any deviations or adaptations made in the use of the standard?	
Were deviations or adaptations made beyond what is specified in the FDA SIS?	×
Were there any exclusions from the standard?	X
Is there an FDA guidance ^a that is associated with this standard?	
 ¹ The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication] ² Authority [21 U.S.C. 360d], http://www fda.gov/MedicalDevices/ DeviceRegulationandGuidence/Standards/default.htm ³ Authority [21 U.S.C. 360d], http://www fda.gov/MedicalDevices/ DeviceRegulationandGuidence/Standards/default.htm ⁴ The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of mathods are described, deviations from the 	tion on

i,

Page 1

EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE						
STANDARD TITLE ISO 10993-10:2002/A	STANDARD TITLE ISO 10993-10:2002/Amd. 1:2006(E) Biological Evaluation of Medical Devices - Part 10: Tests for Irritation and Delayed-Type					
	CONFORMANCE	WITH STANDARD SECTIONS*				
SECTION NUMBER All Sections	SECTION TITLE Various		CONFORMANCE?	/A		
	R OPTION SELECTED * include the guinea pig maximization t	est				
DESCRIPTION						
JUSTIFICATION						
SECTION NUMBER	SECTION TITLE					
TYPE OF DEVIATION O	R OPTION SELECTED •					
DESCRIPTION						
JUSTIFICATION						
SECTION NUMBER	SECTION TITLE		CONFORMANCE?			
TYPE OF DEVIATION O	R OPTION SELECTED *			γ <u>Α</u>		
DESCRIPTION			,			
JUSTIFICATION	<u>_</u>	·····				
explanation is neede described and adequ selected when follow report. More than or * Types of deviations a	d under "justification." Some standard pately justified as appropriate for the so ring a standard is required under "type he page may be necessary. can include an exclusion of a section i	ate whether conformance is met. If a section is include options, so similar to deviations, to ubject device. Explanation of all deviations of deviation or option selected," "description in the standard, a deviation brought out by to to the device, or any adaptation of a section	the option chosen needs to b or description of options on" and "justification" on the the FDA supplemental	be		
	Paperwork I	Reduction Act Statement				
time for review completing and	ing instructions, searching existing da	ation is estimated to average 1 hour per res ata sources, gathering and maintaining the ion. Send comments regarding this burden ggestions for reducing this burden to:	data needed, and			
Food a Office 1350 1	tment of Health and Human Services and Drug Administration of Chief Information Officer Piccard Drive, Room 400 ville, MD 20850	An agency may not conduct or spo required to respond to, a collectio, displays a currently valid OMB co	n of information unless it			
ORM EDA 3654 (6/11		Page 2				

Form Approved: OMB No. 0910-0120; Expiration Date: 12/31/13					
Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)					
This report and the Summary Report Table are to be completed by the applicant when submitting a ences a national or international standard. A separate report is required for each standard referenced i	510(k) ti n the 51	hat refer- 0(k).			
TYPE OF 510(K) SUBMISSION					
STANDARD TITLE ' ANSI/AAMI/ISO 10993-11:2006 Biological Evaluation of Medical Devices - Part 11: Tests for Systemic Toxicity	1				
Please answer the following questions	Yes	No			
Is this standard recognized by FDA ² ?	X				
FDA Recognition number ³	<u></u> 2-118				
Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)?	8				
Is a summary report ⁴ describing the extent of conformance of the standard used included in the 510(k)? If no, complete a summary report table.		X			
Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device?	×				
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).	X				
Does this standard include more than one option or selection of tests? If yes, report options selected in the summary report table.	×				
Were there any deviations or adaptations made in the use of the standard? If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) ⁵ ?					
Were deviations or adaptations made beyond what is specified in the FDA SIS? If yes, report these deviations or adaptations in the summary report table.		X			
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.		X			
Is there an FDA guidance ⁶ that is associated with this standard? If yes, was the guidance document followed in preparation of this 510k? Title of guidance:					
 ¹ The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication] ² Authority [21 U.S.C. 360d], http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm ³ http://www.accessdata.fda.gov/scripts/odm/cfdocs/cfStandards/search.cfm ⁴ The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard: requirements not applicable to the device, and the name and ⁴ The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard: requirements not applicable to the device, and the name and 	udes information by ice. I information und at http inds/search n be found	mation on on which // n c fm Jat			

1

Page 1

EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE						
STANDARD TITLE ANSI/AAMI/ISO 10993-11:2006 Biological Evaluation of Medical Devices - Part 11: Tests for Systemic Toxicity						
	CONFORMANCE W	ITH STANDARD SECTIONS*				
SECTION NUMBER	SECTION TITLE		CONFORMANCE?			
All Sections	Various		🔀 Yes 🗌 No	[]] N/A		
		ty study in rats following subcutaneous im	plantation and sys	temic		
JUSTIFICATION						
SECTION NUMBER	SECTION TITLE		CONFORMANCE?	,		
			Yes No	□ N/A		
TYPE OF DEVIATION	OR OPTION SELECTED *		1			
DESCRIPTION						
JUSTIFICATION			<u>.</u>			
SECTION NUMBER	SECTION TITLE		CONFORMANCE?			
			Yes No	N/A		
TYPE OF DEVIATION	DR OPTION SELECTED +					
DESCRIPTION						
JUSTIFICATION	······································		_			
explanation is need described and adec selected when follo report. More than c * Types of deviations	ed under "justification." Some standards juately justified as appropriate for the sul wing a standard is required under "type o one page may be necessary. I can include an exclusion of a section in	te whether conformance is met. If a section include options, so similar to deviations, th oject device. Explanation of all deviations of of deviation or option selected," "description the standard, a deviation brought out by th o the device, or any adaptation of a section.	e option chosen no r description of op n" and "justification le FDA supplemen	eeds to be tions " on the		
	Paperwork R	eduction Act Statement				
time for revie completing an	ng burden for this collection of informat wing instructions, searching existing dat	ion is estimated to average 1 hour per resp a sources, gathering and maintaining the d n. Send comments regarding this burden e	ata needed, and			
Food Offic 1350	artment of Health and Human Services and Drug Administration se of Chief Information Officer Piccard Drive, Room 400 sville, MD 20850	An agency may not conduct or spon required to respond to, a collection displays a currently valid OMB con	of information unles			
	<u> </u>					

Form Approved: OMB No. 0910-0120, Expiration Date: 12/31/13						
Food and Drug STANDARDS DATA	n and Human Services 9 Administration REPORT FOR 510(k)s n by applicant)					
This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).						
TYPE OF 510(K) SUBMISSION						
STANDARD TITLE ' ANSI/AAMI/ISO 11135-1: 2007 Sterilization of Health Care Prod	ucts – Ethylene Oxide Part 1: Requirements fo	r Develo	pment			
Please answer the following questions		Yes	No			
Is this standard recognized by FDA 2?		×.				
FDA Recognition number ^a		#14-228				
Was a third party laboratory responsible for testing conform in the 510(k)?			X			
Is a summary report ⁴ describing the extent of conformance 510(k)? If no, complete a summary report table			\boxtimes			
Does the test data for this device demonstrate conformity to pertains to this device?		X				
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).		X				
Does this standard include more than one option or selectio If yes, report options selected in the summary report table.	n of tests?		$\overline{\mathbf{X}}$			
Were there any deviations or adaptations made in the use of If yes, were deviations in accordance with the FDA supplement						
Were deviations or adaptations made beyond what is specified of the specified of the section of the summar of the summar sections in the summar sections of the summar sections of the summar sections of the section of			X			
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			X			
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation of Title of guidance:	of this 510k?					
 The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication] Authority [21 U S.C. 360d]. http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm 	address of the test laboratory or certification body invi- assessment to this standard. The summary report inc all standards utilized during the development of the di The supplemental information sheet (SIS) is additional	ludes infor evice il informati	mation on on which			
 ³ http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm ⁴ The summary report should include, any adaptations used to adapt to the device under review (for example, alternative test methods), choices made when options or a selection of methods are described, deviations from the standard; requirements not applicable to the device, and the name and 	is necessary before FDA recognizes the standard. Fo www accessdata.fda gov/scnpts/cdrh/c/docs/cfStanda , The online search for CDRH Guidance Documents ca http://www.fda.gov/MedicalDevices/DeviceRegulation GuidanceDocuments/default.htm	ards/search	n cfm tat			

OSC Phylodiane Service of the Observer, 13

5

		NDARD CONFORMANCE	· ·		
STANDARD TITLE ANSI/AAMI/ISO T	1135-1: 2007 Sterilization of Health Care	Products – Ethylene Oxide Part 1: Requirer	ments for	Develop	ment
	CONFORMANCE W	ITH STANDARD SECTIONS*			
SECTION NUMBER	SECTION TITLE Various		CONFORM	ANCE?	
	OR OPTION SELECTED *	L			
DÉSCRIPTION					
JUSTIFICATION					
SECTION NUMBER	SECTION TITLE	[CONFORM	ANCE?	11 TT
TYPE OF DEVIATION	OR OPTION SELECTED *		Yes	No	<u>[]</u> N/A
DESCRIPTION					
JUSTIFICATION					
2507.000 Million					
SECTION NUMBER	SECTION TITLE				□ N/A
TYPE OF DEVIATION	OR OPTION SELECTED *				
DESCRIPTION					
JUSTIFICATION					· · · · · · · · ·
explanation is need described and ade selected when follo report. More than * Types of deviations	ted under "justification." Some standards quately justified as appropriate for the sub wing a standard is required under "type o one page may be necessary. s can include an exclusion of a section in	te whether conformance is met. If a section include options, so similar to deviations, the oject device. Explanation of all deviations or of deviation or option selected," "description" the standard, a deviation brought out by the the device, or any adaptation of a section	option chi description and "justif	osen nee n of optic ication" o	eds to be ons on the
	Paperwork R	eduction Act Statement			
time for revie completing a	wing instructions, searching existing data	ion is estimated to average 1 hour per resport a sources, gathering and maintaining the data n. Send comments regarding this burden esti- gestions for reducing this burden to:	a needed.	and	
Fou Offic 1350	artment of Health and Human Services d and Drug Administration ce of Chief Information Officer (Piccard Drive, Room 400 kville, MD 20850	An agency may not conduct or sponso required to respond to, a collection of displays a currently valid OMB contro	t informatie		

Form Approved_OMB No. 0910-0120, Expiration Date: 12/31/13					
Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)					
This report and the Summary Report Table are to be comp ences a national or international standard. A separate report					
TYPE OF 510(K) SUBMISSION	Abbreviated				
STANDARD TITLE ANSI/AAMI/ISO 11137-1: 2006 Sterilization of health care produ	ucts – Radiation – Requirements for developme	ent, valio	lation		
Please answer the following questions		Yes	No		
Is this standard recognized by FDA 2?		X			
FDA Recognition number ³		±14-297			
Was a third party laboratory responsible for testing conformi in the 510(k)?			X		
Is a summary report 4 describing the extent of conformance 510(k)? If no, complete a summary report table.			X		
Does the test data for this device demonstrate conformity to pertains to this device?		×			
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).		×			
Does this standard include more than one option or selection If yes, report options selected in the summary report table.	n of tests?		X		
Were there any deviations or adaptations made in the use o If yes, were deviations in accordance with the FDA supplem					
Were deviations or adaptations made beyond what is specified of the summary of the summary set of the set of the summary set of the set of			X		
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			X		
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation o Title of guidance:	of this 510k?		×		
 ¹ The formatting convention for the title is [SDO] [numeric identifier] (title of standard] [date of publication] ² Authonty [21 U.S.C. 360d], http://www.ida.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm ³ http://www.accessdata.fda.gov/scnpts/cdrh/cfdocs/cfStandards/search.cfm ⁴ The summary report should include: any adaptations used to adapt to the 	address of the test laboratory or certification body invo assessment to this standard. The summary report incl all standards utilized during the development of the de * The supplemental information sheet (SIS) is additional is necessary before FDA recognizes the standard. For www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStanda The online search for CDRH Guidance Documents ca	udes infon avice. I informatio und at http rds/search	mation on on which of/ n cfm		
device under review (for example, alternative lest methods), choices made when options or a selection of methods are described; deviations from the standard, requirements not applicable to the device, and the name and	The online search for CDRH Guidance Documents ca http://www.fda.gov/MedicalDevices/DeviceRegulations GuidanceDocuments/default htm	andGuidar	nce/		

IN Rolating Server (open Green TY

		ANDARD CONFORMANCE RY REPORT TABLE			
STANDARD TITLE • ANSI/AAMI/ISO 11	137-1: 2006 Sterilization of health ca	re products – Radiation – Requirements for	developme	nt. valida	ation
	CONFORMANCE	WITH STANDARD SECTIONS*	<u> </u>		
SECTION NUMBER	SECTION TITLE		CONFORM	IANCE?	
All	Various		Yes Yes	No No	
	R OPTION SELECTED * all indicated sections.				
DESCRIPTION	<u> </u>			- N.	
JUSTIFICATION	<u></u>				
SECTION NUMBER	SECTION TITLE	· · · · · · · · · · · · · · · · · ·	CONFORM	ANCE?	
			📋 Yes	No No	🗌 N/A
TYPE OF DEVIATION O	R OPTION SELECTED *		200	51 51. 	41
DESCRIPTION					
JUSTIFICATION					
SECTION NUMBER	SECTION TITLE		CONFORM	ANCE?	
			Yes	No No	□ N/A
TYPE OF DEVIATION O	R OPTION SELECTED *				
DESCRIPTION					
JUSTIFICATION		···- ·····			
explanation is neede described and adequ selected when follow report. More than or * Types of deviations	ed under "justification." Some standard pately justified as appropriate for the se ving a standard is required under "type he page may be necessary. can include an exclusion of a section i	ate whether conformance is met. If a section s include options, so similar to deviations, th ubject device. Explanation of all deviations of of deviation or option selected," "description in the standard, a deviation brought out by th	e option ch or descriptio n° and "justi ne FDA supp	osen nee n of optic fication" (eds to be ons on the
Information sneet (S	IS), a deviation to adapt the standard t	to the device, or any adaptation of a section.			
		Reduction Act Statement	_		
time for review completing and	ing instructions, searching existing de	ation is estimated to average 1 hour per resp at sources, gathering and maintaining the d on. Send comments regarding this burden e ggestions for reducing this burden to:	ata needed.	and	
Food : Office 1350	tment of Health and Human Services and Drug Administration of Chief Information Officer Piccard Drive, Room 400 fille, MD 20850	An agency may not conduct or spor required to respond to, a collection displays a currently valid OMB con	of information	on unless	
FORM FDA 3654 (6/11		Page 2		-	

Form Approved. OMB No. 0910-0120; Expiration Date: 12/31/13

Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)				
This report and the Summary Report Table are to be comp ences a national or international standard. A separate report				
TYPE OF 510(K) SUBMISSION				
X Traditional Special	Abbreviated			
STANDARD TITLE 1 ANSI/AAMI/ISO 11137-2: 2009 Sterilization of health care produ	icts - Radiation - Establishing the sterilization	dose		
Please answer the following questions		Yes	No	
Is this standard recognized by FDA 2?		\boxtimes		
FDA Recognition number ³		# <u>14-225</u>		
Was a third party laboratory responsible for testing conform in the 510(k)?			×	
Is a summary report ⁴ describing the extent of conformance 610(k)? If no, complete a summary report table.				
Does the test data for this device demonstrate conformity to pertains to this device?		×		
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).		×		
Does this standard include more than one option or selection of selection of the summary report table.		×		
Were there any deviations or adaptations made in the use of If yes, were deviations in accordance with the FDA suppler				
Were deviations or adaptations made beyond what is speci If yes, report these deviations or adaptations in the summar			×	
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.				
Is there an FDA guidance ⁶ that is associated with this stan If yes, was the guidance document followed in preparation Title of guidance:	of this 510k?			
 The formatting convention for the title is. [SDO] [numeric identifier] [title of standard] [date of publication] 	address of the test laboratory or certification body in assessment to this standard. The summary report in all standards utilized during the development of the	ncludes inform	normance nation on	
 Authority [21 U.S.C. 360d]. http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm all standards utilized during the development of the device provide the supplementation sheet (SIS) is additional 			on which	
http://www.accessdata.fda.gov/scnpts/cdrh/cfdocs/cfStandards/search cfm	is necessary before FDA recognizes the standard. I www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStan	dards/search	ı cîm	
* The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and	8 The online search for CDRH Guidance Documents http://www.Ida.gov/MedicalDevices/DeviceRegulatin GuidanceDocuments/default.htm			
EORM EDA 3654 (6/11) Par		-		

RM FDA 3654 (6/11)

٢

:

EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE						
STANDARD TITLE ANSI/AAMI/ISO 111	STANDARD TITLE ANSI/AAMI/ISO 11137-2: 2009 Sterilization of health care products - Radiation - Establishing the sterilization dose					
	CONFORMANCE WI	TH STANDARD SECTIONS"		<u></u>		
SECTION NUMBER	SECTION TITLE		CONFORM	ANCE?		
1-6, 9-11	Various		🔀 Yes	No []		
	R OPTION SELECTED • all indicated sections.		L			
DESCRIPTION						
JUSTIFICATION						
SECTION NUMBER	SECTION TITLE		CONFORM	ANCE?		
7-8	Dose Setting using Bioburden/Fraction	Positive Information	☐ Yes	No No	🔀 N/A	
TYPE OF DEVIATION O	R OPTION SELECTED *		<u></u>			
DESCRIPTION			<u></u>			
JUSTIFICATION Methods described in	sections 7 and 8 are not utilized; VDmax	method is utilized as described in sectio	n 9.			
SECTION NUMBER	SECTION TITLE		CONFORM	IANCE?		
	1		🗌 Yes	门 No	🗌 N/A	
TYPE OF DEVIATION O	R OPTION SELECTED *					
DESCRIPTION						
JUSTIFICATION				<u>-</u>		
explanation is neede described and adequ selected when follow report. More than or * Types of deviations of	d under "justification." Some standards in lately justified as appropriate for the subjuing a standard is required under "type of the page may be necessary.	whether conformance is met. If a section include options, so similar to deviations, th ect device. Explanation of all deviations o deviation or option selected," "description ne standard, a deviation brought out by th he device, or any adaptation of a section	e option ch r descriptio 1° and "justi	osen nee n of optic fication" (ads to be ons on the	
	Paperwork Rec	luction Act Statement		===		
time for review completing and	g burden for this collection of informatio ing instructions, searching existing data	n is estimated to average 1 hour per resp sources, gathering and maintaining the da . Send comments regarding this burden e	ata needed,	and		
Food a Office 1350 I	tment of Health and Human Services and Drug Administration of Chief Information Officer Piccard Drive, Room 400 fille, MD 20850	An agency may not conduct or spon required to respond to, a collection displays a currently valid OMB con	of informati	on unless		

	Form Approved: OMB No. 0910-0120; Exp	piration Da	ate: 12/31/13
Food and Drug STANDARDS DATA	and Human Services g Administration REPORT FOR 510(k)s n by applicant)		
This report and the Summary Report Table are to be comp ences a national or international standard. A separate report			
TYPE OF 510(K) SUBMISSION		<u> </u>	
🔀 Traditional 🔄 Special	Abbreviated		
STANDARD TITLE · · ASTM F136-08 Standard Specification for Wrought Titanium-6/	Aluminum-4Vanadium ELI (Extra Low Interstit	ial) Allo	y for
Please answer the following questions		Yes	No
Is this standard recognized by FDA ² ?		X	
FDA Recognition number ³		¥8-164	
Was a third party laboratory responsible for testing conform in the 510(k)?			X
Is a summary report 4 describing the extent of conformance 510(k)? If no, complete a summary report table.			X
Does the test data for this device demonstrate conformity to pertains to this device?		X	
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).		×	
Does this standard include more than one option or selection If yes, report options selected in the summary report table.	n of tests?		X
Were there any deviations or adaptations made in the use of If yes, were deviations in accordance with the FDA supplem			
Were deviations or adaptations made beyond what is speci- If yes, report these deviations or adaptations in the summar			\mathbf{X}
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.		-	X
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation of Title of guidance:	of this 510k?		
 The formatting convention for the title is [SDO] [numeric identifier] [title of standard] [date of publication] Authonty [21 U.S.C. 360d], http://www.fda.gov/MedicalDevices/ 	address of the test laboratory or certification body invo assessment to this standard. The summary report incl all standards utilized during the development of the de	udes infor	nformance mation on
DeviceRegulationandGuidance/Standards/default htm	5 The supplemental information sheet (SIS) is additional is necessary before FDA recognizes the standard. For		
³ http://www.accessdata.tda.gov/scripts/cdm/cfdocs/cfStandards/search.cfm 4.The superson second should usefulde: apti adaptations used to edapt to the second second	www.accessdata.fda.gov/scnpts/cdrh/cfdocs/cfStanda	rds/search	n.cfm
The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods), choices made when options or a selection of methods are described; deviations from the standard, requirements not applicable to the device; and the name and	The online search for CDRH Guidance Documents ca http://www.fda.gov/MedicalDevices/DeviceRegulationa GuidanceDocuments/default.htm	n be found andGuidar	d at nce/

Ind metabolic herodas (Mirakovani) - PR

[NDARD CONFORMANCE			
STANDARD TITLE ASTM F136-08 Sta	undard Specification for Wrought Titanic	m-6Aluminum-4Vanadium ELI (Extra Lo	w Interstitial) Alloy	for
	CONFORMANCE W	ITH STANDARD SECTIONS*			
SECTION NUMBER	SECTION TITLE Various		CONFORMA	NCE?	1 N/A
TYPE OF DEVIATION OR OPTION SELECTED * None. Conformance to all indicated sections.					
DESCRIPTION					
JUSTIFICATION					
SECTION NUMBER	SECTION TITLE			NCE? ☐ Nø	[_] N/A
TYPE OF DEVIATION	OR OPTION SELECTED +				
DESCRIPTION				и .	
JUSTIFICATION					
SECTION NUMBER	SECTION TITLE		CONFORMA		/~~
TYPE OF DEVIATION (DR OPTION SELECTED *		Yes [_] No	N/A
DESCRIPTION					
JUSTIFICATION					
* For completeness list all sections of the standard and indicate whether conformance is met. If a section is not applicable (N/A) an explanation is needed under "justification." Some standards include options, so similar to deviations, the option chosen needs to be described and adequately justified as appropriate for the subject device. Explanation of all deviations or description of options selected when following a standard is required under "type of deviation or option selected," "description" and "justification" on the report. More than one page may be necessary.					
* Types of deviations can include an exclusion of a section in the standard, a deviation brought out by the FDA supplemental information sheet (SIS), a deviation to adapt the standard to the device, or any adaptation of a section.					
Paperwork Reduction Act Statement					
Public reporting burden for this collection of information is estimated to average 4 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to:					
Food Offic 1350	irtment of Health and Human Services and Drug Administration the of Chief Information Officer Piccard Drive, Room 400 wille, MD 20850	An agency may not conduct or spon- required to respond to, a collection displays a currently valid OMB com	of information		

Ĩ.

Form Approved, OMB No. 0910-0120; Expiration Date: 12/31/13	Form Approved.	OMB No.	0910-0120:	Expiration Date:	12/31/13
---	----------------	---------	------------	------------------	----------

Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)					
This report and the Summary Report Table are to be com ences a national or international standard. A separate report	pleted by the applicant when submitting a rt is required for each standard referenced	510(k) t in the 5	hat refer- 10(k).		
TYPE OF 510(K) SUBMISSION					
🔀 Traditional 🔄 Special	Abbreviated				
STANDARD TITLE ' ASTM F648-07 Standard Specification for Ultra-High-Molecular-	Weight Polyethylene Powder and Fabricated Fe	orm for S	Surgica)		
Please answer the following questions		Yes	No		
Is this standard recognized by FDA 2?		×			
FDA Recognition number ³		# 8-178			
Was a third party laboratory responsible for testing conform in the 510(k)?	ity of the device to this standard identified		×		
Is a summary report ⁴ describing the extent of conformance 510(k)? If no, complete a summary report table.			×		
Does the test data for this device demonstrate conformity to pertains to this device?		X			
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).		X			
Does this standard include more than one option or selection If yes, report options selected in the summary report table.	on of tests?				
Were there any deviations or adaptations made in the use of If yes, were deviations in accordance with the FDA supplement					
Were deviations or adaptations made beyond what is speci If yes, report these deviations or adaptations in the summar			X		
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			×		
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation Title of guidance:	of this 510k?				
 ¹ The formatting convention for the title is [SDO] [numeric identifier] [title of standard] (date of publication] ² Authority [21 U.S.C. 360d], http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm ³ http://www.accessdata.fda.gov/scnpts/cdrh/cfdocs/cfStandards/search.cfm 	address of the test laboratory or certification body inva assessment to this standard. The summary report inc all standards utilized during the development of the de s The supplemental information sheet (SIS) is additional is necessary before FDA recognizes the standard. Fo www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStanda	ludes infor evice al informati und at http	mation on on which off		
* The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described, deviations from the standard, requirements not applicable to the device; and the name and	6 The online search for CDRH Guidance Documents or http://www.fda.gov/MedicalDevices/DeviceRegulation GuidanceDocuments/default.htm				

1

Page 1

EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE				
STANDARD TITLE ASTM F648-07 Stands	ard Specification for Ultra-High-Molecular-We	ight Polyethylene Powder and Fal	pricated Form for Surgical	
	CONFORMANCE WITH STA	NDARD SECTIONS*		
SECTION NUMBER	SECTION TITLE		CONFORMANCE?	
All	Various		🕅 Yes 🗌 No 🗍 N/A	
TYPE OF DEVIATION OF None. Conformance to	all indicated sections.			
DESCRIPTION				
JUSTIFICATION				
SECTION NUMBER	SECTION TITLE		CONFORMANCE?	
			Yes No N/A	
DESCRIPTION				
JUSTIFICATION				
SECTION NUMBER	SECTION TITLE		CONFORMANCE?	
			Yes No N/A	
TYPE OF DEVIATION OF	ROPTION SELECTED *			
DESCRIPTION		· · · · · · · · · · · · · · · · · · ·		
JUSTIFICATION		_		
		_		
explanation is needed described and adequ selected when follow	all sections of the standard and indicate whether d under "justification." Some standards include of ately justified as appropriate for the subject dev ing a standard is required under "type of deviation e page may be necessary.	options, so similar to deviations, th ice. Explanation of all deviations o	e option chosen needs to be r description of options	
	an include an exclusion of a section in the stand S), a deviation to adapt the standard to the devi		e FDA supplemental	
	Paperwork Reduction	Act Statement	فاعتلدته وينبغه فبمتعد والمتعادي	
time for review completing and	g burden for this collection of information is est ing instructions, searching existing data sources reviewing the collection of information. Send ollection of information, including suggestions	imated to average 1 hour per resp s, gathering and maintaining the d comments regarding this burden e	ata needed, and	
Food a Office 1350 F	ment of Health and Human Services and Drug Administration of Chief Information Officer Piccard Drive, Room 400 ille, MD 20850	An agency may not conduct or spon required to respond to, a collection displays a currently valid OMB con	of information unless it	

	Form Approved: OMB No. 0910-0120: Ex	piration Da	ate. 12/31/13
Food and Drug STANDARDS DATA	n and Human Services g Administration REPORT FOR 510(k)s n by applicant)		
This report and the Summary Report Table are to be comp ences a national or international standard. A separate report			
TYPE OF 510(K) SUBMISSION			
ASTM F1537-00 Specification for Wrought Cobalt-28 Chromiur	n-6 Molybdenum Alloy for Surgical Implants		
Please answer the following questions		Yes	No
Is this standard recognized by FDA 2?		X	
FDA Recognition number ³		#8-182	
Was a third party laboratory responsible for testing conform in the 510(k)?			×
Is a summary report ⁴ describing the extent of conformance 510(k)? If no, complete a summary report table.			X
Does the test data for this device demonstrate conformity to pertains to this device?		×	
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).		X)	
Does this standard include more than one option or selection of selection of selected in the summary report table	on of tests?		X
Were there any deviations or adaptations made in the use of If yes, were deviations in accordance with the FDA supplen			
Were deviations or adaptations made beyond what is speci If yes, report these deviations or adaptations in the summar			X
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation of Title of guidance:	of this 510k?		
 ¹ The formatting convention for the title is [SDO] (numeric identifier) [title of standard] (date of publication] ² Authority [21 U S C 360d], http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm ³ http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm 	address of the test laboratory or certification body inv assessment to this standard. The summary report inc all standards utilized during the development of the di 5 The supplemental information sheet (SIS) is additional is necessary before FDA recognizes the standard. Fo www.accessdata.tda.gov/scnpts/cdm/c/docs/cfStanda	ludes infor evice. al information of http: ound at http:	mation on on which or//
* The summary report should include any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and	c The online search for CDRH Guidance Documents ca http://www.fda.gov/MedicalDevices/DeviceRegulation GuidanceDocuments/default htm	an be found	d al

.

1

PSC Publishing Services (101), 4147 (201), FF

1

EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE				
STANDARD TITLE • ASTM F1537-00 S	pecification for Wrought Cobalt-28 Ch	romium-6 Molybdenum Alloy for Surgical Implants		
	CONFORMANCE	WITH STANDARD SECTIONS"		
SECTION NUMBER	SECTION TITLE Various			
TYPE OF DEVIATION None. Conformance	OR OPTION SELECTED *			
DESCRIPTION				
JUSTIFICATION				
SECTION NUMBER	SECTION TITLE	CONFORMANCE?		
TYPE OF DEVIATION	OR OPTION SELECTED .			
DESCRIPTION				
JUSTIFICATION				
SECTION NUMBER	SECTION TITLE			
TYPE OF DEVIATION	OR OPTION SELECTED +			
DESCRIPTION				
JUSTIFICATION				
explanation is need described and aded selected when follo report. More than o * Types of deviations	led under "justification." Some standard quately justified as appropriate for the si wing a standard is required under "type one page may be necessary. s can include an exclusion of a section i	ate whether conformance is met. If a section is not applicable (N/A) an s include options, so similar to deviations, the option chosen needs to be ubject device. Explanation of all deviations or description of options of deviation or option selected," "description" and "justification" on the in the standard, a deviation brought out by the FDA supplemental		
information sheet ()		to the device, or any adaptation of a section.		
time for revie completing ar	ng burden for this collection of information wing instructions, searching existing data	Reduction Act Statement ation is estimated to average 1 hour per response, including the ata sources, gathering and maintaining the data needed, and on. Send comments regarding this burden estimate or any other agestions for reducing this burden to:		
Depa Fuoc Offic 1350	artment of Health and Human Services I and Drug Administration Se of Chief Information Officer Priceard Drive, Room 400 Wille, MD 20850	An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number		



Device Name

Zimmer Nexel Total Elbow

Reason for the Submission

Zimmer is seeking clearance to market the Zimmer Nexel Total Elbow, which has not previously been marketed.

Device Description

Overview

Constant of relation The Zimmer Nexel Total Elbow is designed to replace the articulating surfaces of and restore motion to the humeroulnar (elbow) joint. The implant is a semi-constrained device assembly and consists of a humeral component, an ulnar component, a---humeral bearing-A, 2 ulnar bearings-B, an axle pin and 2 humeral screws.




Image 10.1: The Zimmer Nexel Total Elbow, exploded view

A detailed description of the proposed device is presented in Section 11, Device Description.



Indications for Use

Indications for use include:

- Elbow joint destruction which significantly compromises the activities of daily living
- Post-traumatic lesions or bone loss contributing to elbow instability
- Ankylosed joints, especially in cases of bilateral ankylosis from causes other than active sepsis
- Advanced rheumatoid, post-traumatic, or degenerative arthritis with incapacitating pain
- Instability or loss of motion when the degree of joint or soft tissue damage precludes reliable osteosynthesis
- Acute comminuted articular fracture of the elbow joint surfaces that precludes less radical procedures, including 13-C3 fractures of the distal humerus
- Revision arthroplasty

CAUTION: This device is intended for cemented use only.

Submission Overview

The table below provides a comparison between the proposed and predicate devices. Additional information is presented in <u>Section 12</u>, <u>Substantial Equivalence</u> <u>Discussion</u>.

		Proposed Device: The Zimmer Nexel Total Elbow	Predicate Device: The Coonrad/Morrey Total Elbow (K001918, K053189)
	Stem with Triangular Cross-Section	X	x
T E	Straight Stem	x	x
Humeral Component	Anterior Flange	X	x
and m	Titanium Alloy	X	x
Ŭ	Titanium Beads on Distal Region		x
	Titanium Plasma Spray on Distal Region	X	
1	Stem with Quadrangular Cross-Section	X	x
ent	Curved Stem	X	X
Ulnar Component	Nitrogen Enriched Surface Hardening of Ulnar Eye	x	
S S	Titanium Alloy	X	X
	Titanium Plasma Spray on Proximal Region	x	x
20	2-Part, Hollow Snap-Pin		X
Linkage	Solid Axle Pin	X	
- ご	Screws	X	
	Semi-Constrained	X	X
2.5	Terminal Sterilization	X	X
	Titanium-on-Polyethylene Articulation	<u>x</u>	X



Performance Testing

The following performance testing was conducted on the Zimmer Nexel Total Elbow. A detailed summary and copy of each test report are included in <u>Section 18, Performance</u> <u>Testing – Bench</u>.

Stem Fatigue Testing

The fatigue strength of the Zimmer Nexel Total Elbow humeral component distal stem at (b)(4) load cycles is (b)(1, which is (b) stronger than that of the legally marketed Zimmer Coonrad/Morrey Total Elbow humeral component (b)(4) at a comparable point on the distal stem.

The fatigue strength of the Zimmer Nexel Total Elbow ulnar component stem in the mid-stem region at (b)(4) load cycles is at least (b)(c) (the fatigue strength of the legally marketed Zimmer Coonrad/Morrey Total Elbow ulnar component at a comparable point on the stem).

The fatigue strength of the Zimmer Nexel Total Elbow ulnar component stem in the plasma spray (proximal stem) region at (b)(4) load cycles is at least (b)(c) (the fatigue strength of the legally marketed Zimmer Coonrad/Morrey Total Elbow ulnar component at a comparable point on the stem).

These studies demonstrated that the Zimmer Nexel Total Elbow is at least as safe and effective as the predicate device (Zimmer Coonrad/Morrey Total Elbow) in terms of humeral and ulnar stem fatigue strength.

Wear Testing

The Zimmer Nexel Total Elbow mean gravimetric wear rate (b)(4) at (b)(4) at (b)(4) load cycles) was (b) less than that of the predicate.

This study demonstrated that the Zimmer Nexel Total Elbow met or exceeded all acceptance criteria and is at least as safe and effective as the predicate device in terms of wear performance.

Durability Testing

The Zimmer Nexel Total Elbow device achieved (b)(4) prun-outs at the equivalent of (b)(4) weight-in-hand(b)(4) peak JRF and(b)(4) varus-valgus moment), a load level (b)(higher than the load level at which the predicate device ran out.

This study demonstrated that the Zimmer Nexel Total Elbow is at least as safe and effective as the predicate device in terms of durability.



Modular Connection Fatigue Testing

The median fatigue strength of the non-articulating, mechanically locked, modular Zimmer Nexel Total Elbow implant components was no less than (b)(4)

This study demonstrated that the Zimmer Nexel Total Elbow is safe and effective in terms of non-articulating, mechanically locked, modular connection fatigue.

Support Research

The following internal (Zimmer) research and testing informed and/or support the aforementioned testing and can be found in their entirety at the end of <u>Section 18</u>, <u>Performance Testing – Bench</u>.

- Zimmer Research Memo ZRM_WA_0179_09, Analysis of Retrieved <u>Coonrad/Morrey Elbow Total Arthroplasty Components</u>
- Zimmer Research Memo ZRM_WA_0193_09, Development of the Coonrad/Morrey Elbow Durability Test
- Zimmer Research Report ZRR_WA_1988_09, Literature Review Summary of Elbow Joint Mechanics and Total Elbow Arthroplasty Outcomes
- Zimmer Research Report ZRR_WA_2252_10, Evaluation of the AS/BF Taper Adaptor Component Using Electrochemical Test Methods
- Zimmer Technical Memo 1222.00, Fatigue Strength of Plasma Sprayed Ti-6AL-4V Alloy Material
- Zimmer Research Memo ZRM_WA_0292 11, Summary of Clinical Complications of the Coonrad/Morrey Elbow System
- Zimmer Research Report ZRR_WA_2338_11, Vitamin E Highly Crosslinked Polyethylene Shelf Life Sludy
- Zimmer Research Report ZRR_WI_2441_11, Finite Element Analysts of the Coonrad-Morrey Total Elbow System
- Zimmer Research Report ZRR_WA_2476_11, Biomechanics Rationale for the Loading of Total Elbow Prosthesis During In Vitro Simulations
- Zimmer Research Report ZRR_WA_2481_11, Determination of the Functional Torque Range of the New Zimmer Total Elbow Humeral Screw and Hex Driver
- Zimmer Technical Memo 1217.00, Coonrad/Morrey Elbow Ulnar Stem Fatigue Testing
- <u>Zimmer Technical Memo 1218.00</u>, <u>Coonrad/Morrey Humeral Stem Fatigue</u> <u>Testing</u>
- <u>Memo to DHF File Z07-014, Rationale for Determination of Torque</u> <u>Specification and Torque Range Requirement for the Mechanical Fastening</u> <u>Used in Zimmer New Total Elbow Prosthesis</u>



- Zimmer Research Report ZRR_WA_2646_12, Pin-on-Flat Wear Comparison of the New Zimmer Total Elbow and the Coonrad-Morrey Articulation Couples
- Zimmer Research Report ZRR_WA_2648_12; Finite Element Analysis of Nexel Elbow Ulnar Implant Neck Region
- Zimmer Research Report ZRR_WI_2441_11, Finite Element Analysis of the Coonrad-Morrey Total Elbow System



61

Device Description

Overview

The Zimmer Nexel Total Elbow is an implant designed to replace the articulating surfaces of and restore motion to the humeroulnar (elbow) joint. The implant is a semi-constrained device assembly and consists of the following components, which are described in more detail below:

- Humeral Component
- Ulnar Component
- Humeral Bearing-A
- 2 Ulnar Bearings-B
- <u>Axle Pin</u>
- 2 Humeral Screws




Image 11.1: The Zimmer Nexel Total Elbow, exploded view






Image 11.2: The Zimmer Nexel Total Elbow, assembled





Image 11.3: The Zimmer Nexel Total Elbow, assembled with Humeral Component translucent



Humeral Component

The *humeral component* is made from wrought *Tivanium* alloy (Ti-6Al-4V) and has the following major design features: a *stem*, an *anterior flange*, a *yoke* and *plasma spray regions*.



Humeral Component, major design features

The *humeral stem* is intended to be implanted with bone cement into the patient's humeral medullary canal. It has a triangular cross-section, which is designed to prevent rotation of the *humeral component* within the patient's medullary canal; and, it is straight along its proximal-distal axis to conform to the straight shaft of the distal humeral medullary canal.

The *humeral anterior flange* is designed to accept a bone graft, which is wedged between the *anterior flange* and the exterior cortical bone of the humeral shaft. This *anterior flange* and bone graft combination are intended to improve both antirotational stability and posterior migration of the *humeral component*, thereby resisting stem loosening.

The *humeral yoke* has rounded corners (as opposed to angular) on its proximal surface in order to prevent stress risers (that are often associated with bone preparation cuts required for angled corners abutting bone). This rounded *yoke* design is intended to preserve the strength of (and avoid fracture of) the medial and lateral supracondylar columns. The *yoke* walls also house partially threaded holes with Spiralock® technology (tap manufactured by Stanley Black & Decker, Madison Heights, MI), which engage with the *humeral screws*.



The *humeral component* has several regions covered in plasma spray, which is contained within the Plasma Sprayed Titanium-6Al-4V Alloy Coating Material Master File (FDA Ref: MAF-1909, Letter of Authorization can be found in Section 21). The *plasma spray regions* are designed to increase the coefficient of friction between the implant and the bone cement, thereby improving anti-rotational stability and aiding in resistance to movement of the implant with respect to the humerus.

In addition to the aforementioned major design features, the *humeral component* also includes the following features to accept the *axle-pin* and *bearings A and B*: v-grooves at the base of the keyways to accept the *axle-pin*, a hole to accept *humeral bearing-A* and two keyways to accept *ulnar bearings-B*.




Image 11.5: Humeral Component, design features for mating to bearings and axle-pin

Ulnar Component

The *ulnar component* is made from wrought *Tivanium* alloy (Ti-6AI-4V) and has the following major design features: a *stem*, an *eye* and a *plasma spray region*.




Image 11.6: Ulnar Component, major design features



65

The *ulnar stem* is intended to be implanted with bone cement into the patient's ulnar medullary canal. It has a quadrangular cross-section, which is designed to prevent rotation of the *ulnar component* within the patient's medullary canal; and, it is curved along its proximal-distal axis to conform to the curved shaft of the proximal ulnar medullary canal.

The ulnar eye is both highly polished and nitrogen-enriched (as recommended in Section J, Titanium Articular Surface Coatings, of Class II Special Controls Guidance Document: Knee Joint Patellofemorotibial and Femorotibial Metal/Polymer Porous-Coated Uncemented; Guidance for Industry and FDA) to limit wear of the apposing Bearings A and B upon which it articulates.

The proximal portion of the *ulnar stem* is covered with plasma spray, which is contained within the Coatings Characterization of Orthopedic Titanium 6Al4V Coatings Material Master File (FDA Ref: MAF-652, Letter of Authorization can be found in Section 22). The *plasma spray region* is designed to increase the coefficient of friction between the implant and the bone cement, thereby improving anti-rotational stability and aiding in resistance to movement of the *ulnar stem* with respect to the ulna.

The *ulnar component* is side-specific, meaning that there is a left *ulnar component* and a right *ulnar component* available for every size offering.

Humeral Bearing-A

The *Humeral Bearing-A* is made from Vivacit- E^{TM} (Vitamin E [α -tocopherol] stabilized, highly crosslinked ultra-high molecular weight polyethylene), which is contained within the Vitamin E Stabilized HXPE Material Master File authored by Zimmer (FDA Ref: MAF-1868). *Humeral bearing-A* is designed to articulate against the *ulnar component* throughout most of the humeroulnar joint range of motion, and has the following major design features: *concave articulation surface, anterior and posterior ridges*, a *peg*, four *peripheral bumps* and *medial and lateral rails*.






Image 11.7: Humeral Bearing-A, major design features



The humeral bearing-A concave articulation surface is intended to (along with the articulation surface of the ulnar bearings-B) allow for load sharing from the ulnar component onto humeral bearing-A throughout most of the range of motion of the humeroulnar joint. The posterior ridge specifically supports high-flexion motions (see Image 11.8, below). The humeral bearing-A anterior ridge is symmetrical to the posterior ridge in order to remove the human factors issue of having to discern anterior from posterior when assembling humeral bearing-A to the humeral component: humeral bearing-A cannot be put in backwards, because it is symmetrical.





The humeral bearing-A peg is designed to mate to the hole at the proximal base of the humeral yoke with a press (not a snap) fit, which is intended to limit motion between humeral bearing-A and the humeral yoke.

The humeral bearing-A (four) peripheral bumps are designed to mate against the interior walls of the humeral yoke with a press-fit, which is intended to limit motion (and therefore backside wear) of the humeral bearing-A relative to the humeral yoke.

The humeral bearing-A *medial and lateral rails* are designed to allow clearance for (as opposed to a direct mate with) *ulnar bearings-B*, which directly mate to the *humeral yoke* and to the *axle pin*.

Note: The *humeral bearing-A* may also be referred to as "humeral bearing" and "bearing A" throughout this submission.

Ulnar Bearings-B

The *ulnar bearings-B* (two per implant assembly) are made from Vivacit- E^{TM} (Vitamin E [α -tocopherol] stabilized, highly crosslinked ultra-high molecular weight polyethylene), which is contained within the Vitamin E Stabilized HXPE Material

ņ,



Master File authored by Zimmer (FDA Ref: MAF-1868). Ulnar bearings-B are designed to articulate against the *ulnar component* throughout the entire humeroulnar joint range of motion, and have the following major design features: *articulation surface, anti-rotation key, through-hole* and *compressible semi-circle*.




Image 11.9: Ulnar Bearings-B, major design features

The *ulnar bearings-B articulation surface* is intended to (along with the articulation surface of the *ulnar bearing-A*) allow for load sharing from the *ulnar component* onto *humeral bearing-B* throughout the range of motion of the humeroulnar joint.

The *ulnar bearings-B anti-rotation keys* mate into the *humeral yoke keyways* with a press-fit to limit translational and rotational motion of *ulnar bearings-B* with respect to the *humeral component*.

The compressible semi-circle of onc ulnar bearing-B is designed to align with the compressible semi-circle of the second ulnar bearing-B such that a nearly full circle is created. This design feature is intended to ensure proper alignment of ulnar bearings-B with respect to each other (during assembly of ulnar bearings-B to the axle pin using the ulnar bearing assembly tool) and, after proper assembly to the axle pin, to the humeral yoke keyways. Due to the fact that the compressible semi-circle is more easily deformed than the broad, apposing faces of ulnar bearings-B, it creates a press-fit between the ulnar bearings-B and between each ulnar bearing-B and the apposing interior walls of the humeral yoke. This press-fit is intended to limit mediolateral motion of the ulnar bearings-B along the axle-pin, thereby limiting backside wear of the ulnar bearings-B.



The *ulnar bearings-B through-hole* mates around the *axle pin* with a press-fit intended to limit motion of *ulnar bearings-B* with respect to the *axle pin*, and thus limit backside wear of *ulnar bearings-B*. The position of the *through-hole* is eccentric (with respect to the outer diameter of *ulnar bearings-B*) in order to place the largest cross-section of Vivacit-ETM material in opposition to the highest loads (Joint Reaction Forces, JRFs) induced throughout the flexion-extension cycle of the humeroulnar joint (see image 11.10, below).

Note: The *ulnar bearings-B* do NOT mate to or in any way interface with the *humeral screws*.

Note: The *ulnar bearings-B* may also be referred to as "ulnar bearing(s)" and "bearing(s) B" throughout this submission.




Image 11.10: Thickest Vitamin E HXPE placement to oppose highest JRFs

Axle Pin

The *axle pin* is made from *Zimaloy* (Cobalt-Chromium-Molybdenum Alloy), is designed to remain in a fixed position (not translate nor rotate) throughout the entire humeroulnar joint range of motion, and has the following major design feature: *v-groove*.



Image 11.11: Axle Pin, major design feature



69

The *v-groove* is designed to mate with the taper of the *humeral screws* that, upon screw thread engagement in the *humeral yoke*, rigidly compresses the *axle pin* against the *humeral yoke*.

Additionally, the *axle pin* is designed to rigidly mate with *ulnar bearings-B*, the *humeral screws* and the *humeral yoke*, creating a rigid assembly along which joint reaction forces are intended to be shared.

Note: The axle-pin may also be referred to as "axle" throughout this submission.

Humeral Screws

The *humeral screws* (two per implant assembly) are made from *Zimaloy* (Cobalt-Chromium-Molybdenum Alloy), are designed to remain in a fixed position throughout the entire humeroulnar joint range of motion, and have the following major design features: *headless*, *threads*, *pilot* and *taper*.




Image 11.12: Humeral Screws, major design feature

The humeral screw threads engage with threads of the humeral yoke while the humeral screw pilot engages a hole in the wall of the humeral yoke and the humeral screw taper engages with the axle pin v-groove, compressing the ends of the axle pin against the v-groove at the base of the humeral yoke keyway.

The *humeral screws* are *headless*, a feature that allows the *humeral screws* to countersink into the *humeral yoke*, which is intended to ensure the constant application of a clamp-load to the *axle-pin* and to avoid soft tissue impingement or irritation associated with a protruding screw head.

Note: While the *humeral screws* are designed to apply a consistent clamp-load to the *axle pin*, they are *not* directly loaded by any compressive joint reaction force (JRF) during the flexion-extension cycle (this design feature is intended to avoid screw back-out by not allowing any complete transverse load-reversal across the *humeral screw threads*). Instead, compressive JRFs are transmitted from the *ulnar implant*



component to the *ulnar bearings-B* to the *axle pin* to the *v-grooves* at the base of the *humeral yoke keyways* and finally throughout the *humeral component*. See <u>Zimmer</u> <u>Research Report ZRR_WA_2598_12</u>, <u>New Zimmer Total Elbow Prosthesis Modular</u> <u>Connection Fatigue Testing</u> for a more detailed explanation.

System Compatibility

No Zimmer Nexel Total Elbow component is compatible with any device outside the Zimmer Nexel Total Elbow system.

Size Interchangeability

All Zimmer Nexel Total Elbow *humeral component* sizes are designed to mate with all Zimmer Nexel Total Elbow *ulnar component* sizes. The *axle pin* size 4, *humeral bearing-A* size 4, and *ulnar bearings-B* size 4 are designed to mate with all size 4 *humeral components*. The *axle pin* size 5/6, *humeral bearing-A* size 5/6, and *ulnar bearings-B* size 5/6 are designed to mate with all size 5 and size 6 *humeral components*. The *humeral screws* are designed to mate with all Zimmer Nexel Total Elbow *humeral components* and *axle pin* sizes.

Engineering Drawings

:

Representative <u>engineering drawings</u> for the proposed device are included at the end of this section.



Catalog Numbers

All catalog numbers	for the Zimmer	Nexel Total	Elbow are	listed below.
---------------------	----------------	-------------	-----------	---------------

Catalog Number	Description
00-8400-014-07	Zimmer Nexel Total Elbow - Ulnar Stem, Size 4, 75 mm, Left
00-8400-014-11	Zimmer Nexel Total Elbow - Ulnar Stem, Size 4, 115 mm, Left
00-8400-024-07	Zimmer Nexel Total Elbow - Ulnar Stem, Size 4, 75 mm, Right
00-8400-024-11	Zimmer Nexel Total Elbow - Ulnar Stem, Size 4, 115 mm, Right
00-8400-015-07	Zimmer Nexel Total Elbow - Ulnar Stem, Size 5, 75 mm, Left
00-8400-015-11	Zimmer Nexel Total Elbow - Ulnar Stem, Size 5, 115 mm, Left
00-8400-025-07	Zimmer Nexel Total Elbow - Ulnar Stem, Size 5, 75 mm, Right
00-8400-025-11	Zimmer Nexel Total Elbow - Ulnar Stern, Size 5, 115 mm, Right
00-8400-016-09	Zimmer Nexel Total Elbow - Ulnar Stem, Size 6, 90 mm, Left
00-8400-016-11	Zimmer Nexel Total Elbow - Ulnar Stem, Size 6, 115 mm, Left
00-8400-026-09	Zimmer Nexel Total Elbow - Ulnar Stem, Size 6, 90 mm, Right
00-8400-026-11	Zimmer Nexel Total Elbow - Ulnar Stem, Size 6, 115 mm, Right
00-8400-044-10	Zimmer Nexel Total Elbow - Humeral Stem. Size 4, 100 mm
00-8400-044-15	Zimmer Nexel Total Elbow - Humeral Stem, Size 4, 150 mm
00-8400-045-10	Zimmer Nexel Total Elbow - Humeral Stem, Size 5, 100 mm
00-8400-045-15	Zimmer Nexel Total Elbow - Humeral Stem, Size 5, 150 mm
00-8400-055-15	Zimmer Nexel Total Elbow - Humeral Stem, Long Flange, Size 5, 150 mm
00-8400-055-20	Zimmer Nexel Total Elbow - Humeral Stem, Long Flauge, Size 5, 200 mm
00-8400-046-10	Zimmer Nexel Total Elbow - Humeral Stem, Size 6, 100 mm
00-8400-046-15	Zimmer Nexel Total Elbow - Humeral Stem, Size 6, 150 mm
00-8400-056-15	Zimmer Nexel Total Elbow - Humeral Stem, Long Flange, Size 6, 150 mm
00-8400-056-20	Zimmer Nexel Total Elbow - Humeral Stem, Long Flange, Size 6, 200 mm
Kit # 00-8400-094-00	Zimmer Nexel Total Elbow - Articulation Kit, Size 4 Kit Contents:
00-8400-194-01	1 Humeral Bearing-A Size 4
00-8400-194-02 00-8400-194-04	2 Ulnar Bearings- B Size 4 1 Axle Pin Size 4
Kit # 00-8400-095-00	Zimmer Nexel Total Elbow - Articulation Kit, Size 5/6
	Kit Contents:
00-8400-195-01 00-8400-195-02	l Humeral Bearing-A Size 5/6 2 Ulnar Bearings- B Size 5/6
00-8400-195-04	1 Axle Pin Size 5/6
Kit # 00-8400-090-00	Zimmer Nexel Total Elbow - Humeral Screw Kit Kit Contents:
00-8400-190-00	2 Humeral Screws



Materials

The *humeral and ulnar components* are made of wrought *Tivanium* alloy (Ti-6Al-4V) and are titanium plasma sprayed (on regions specified in the <u>engineering drawings</u> at the end of this section). The *bearings A and B* are made of *Vitamin E HXPE* (Vitamin E Highly Crosslinked Polyethylene). The *axle pin* and *humeral screws* are made of *Zimaloy* (Cobalt-Chromium-Molybdenum Alloy).

The following tests contained within the Material Master File for Vitamin E Stabilized Highly Crosslinked Polyethylene (Master File MAF-1868) are applicable to usage in elbow devices and support substantial equivalence of the Vitamin E Stabilized HXPE material when used in elbow devices:

Report Number	Report Name	Master File Location
ZRR_WA_2401_11 Rev, 2	TESTING OF VITAMIN E HIGHLY CROSSLINKED POLYETHYLENE PROPERTIES PER ASTM F2759-09	Amendment 02, Exhibit B
N102082	NAMSA - TOXICOLOGICAL RISK ASSESSMENT	Volume 1, Page 54
T0118_913/S	13 WEEK SYSTEMIC TOXICITY STUDY IN RATS FOLLOWING SUBCUTANEOUS IMPLANTATION	Volume 1, Page 86
V0014_130	CYTOTOXICITY STUDY USING THE ISO ELUTION METHOD	Volume 1, Page 76
V0023_211	BACTERIAL REVERSE MUTATION STUDY - EXTRACTS	Volume 1, Page 108
V0573_000/S	GENOTOXICITY: MOUSE LYMPHOMA ASSAY	Volume 1, Page 125
T0566_500, T0566_501	MOUSE PERIPHERIAL BLOOD MICRONUCLEUS STUDY	Volume 1, Pae 147
T1251_800	ISO INTRACUTANEOUS STUDY IN RABBITS	Volume 1, Page 179
T1261_300	ISO GUINEA PIG MAXIMIZATION SENSITIZATION TEST	Volume 1, Page 189
T1250_812	ISO MUSCLE IMPLANTATION STUDY IN RABBITS - 12 WEEKS	Volume 1, Page 205
T1250_802	ISO MUSCLE IMPLANTATION STUDY IN RABBITS - 2 WEEKS	Volume 1, Page 217
ZRR_WA_2338_11	VITAMIN E HIGHLY CROSSLINKED POLYETHYLENE SHELF LIFE STUDY	Volume 1, Page 230
T02625-500	ISO SYSTEMIC TOXICITY STUDY IN MICE	Volume 1, Page 245
T0118_926	26 WEEK SYSTEMIC TOXICITY STUDY IN RATS FOLLOWING SUBCUTANEOUS IMPLANTATION	Volume 1, Page 256
ZRR_WA_2424_11	EVALUATION OF VITAMIN E DOPED AND CROSSLINKED UHMWPE WEAR DEBRIS IN THE RABBIT KNEE FOLLOWING PERCUTANEOUS INJECTION	Amendment 01, Page 9



73

Report Number	Report Name	Master File Location
10T_60283_02	BIOLOGICAL RESPONSE TO WEAR DEBRIS Note: The particulate study was done to detect local and systemic response to vitamin E grafted UHMWPE. A rabbit model which is widely accepted was used. The particulates were produced from production materials and have size distribution that is representative of those that have been isolated from tissues harvested during hip or knee revision surgeries. The particles were injected directly into the knee capsule of the animals. This study was not focused on any particular implant or device. The emphasis was on inflammatory response to the particulates. The study covered all devices that generate particulates.	Volume 2, Page 243
ZRR_WA_2409_11	VIATMIN E STABILIZED HIGHLY CROSSLINKED POLYETHYLENE EXTENDED AGING STUDY	Volume 2, Page 394
ZRR_WA_2373_11	ENVIRONMENTAL STRESS CRACKING OF VITAMIN E, LONGEVITY, CONVENTIONAL UHMWPE	Volume 2, Page 422
ZRR_WA_2403_11	VITAMIN E STABILIZED HIGHLY CROSSLINKED POLYETHYLENE- EXTRACTABILITY STUDY	Volume 2, Page 441
ZRR_WA_2412_11	VITAMIN E STABILIZED HIGHLY CROSSLINKED - POLYETHYLENE- PRESSURE INDUCED LEACHING	Volume 2, Page 492
ZRR_WA_2382_11	VITAMIN E STABILIZED HIGHLY CROSSLINKED POLYETHYLENE- TRANSFORMATION PRODUCTS EVALUATION BY HPLC AND GC	Volume 2, Page 507
ZRM_WA_0293_11	LAMELLAE THICKNESS	Volume 2, Page 566
ZRM_WA_0277_11	MOLECULAR WEIGHT	Volume 2, Page 578
CPG REPORT_ 11622_1	GC/MS ANALYSIS OF UHMWPE	Amendment 02, Exhibit A

Surface Characteristics

The *eye* of the *ulnar component* is treated with a nitrogen ion implantation (also known as nitriding) process (to increase *Tivanium* alloy surface hardness) on the region specified in the <u>engineering drawings</u> at the end of this section.

The articulating surfaces of the *ulnar component* are polished to an R_a value of 8 microinches or less and the articulating surfaces of the *bearings A and B* are machined to an R_a value of 32 microinches or less to provide for smooth articulating surfaces.

_



74

Mechanism of Action

The Zimmer Nexel Total Elbow implant assembly is a semi-constrained prosthesis designed to allow anatomic reconstruction of the humeroulnar joint with the natural change in humeroulnar joint center of rotation throughout the flexion/extension arc based on research presented in the Zimmer Research Report ZRR_WA_2476_11, Biomechanics Rationale for the Loading of Total Elbow Prosthesis During In Vitro Simulations. See table below for range of motion details.

	Feature				Proposed Device: The Zimmer Nexel Total Elbow
	Varus/Valgus Angular Laxity Internal/External Rotational Laxity				Neutral to Maximum Varus: (b) Neutral to Maximum Valgus: (4) Total Included Angle: (b)
				Neutral to Maximum Internal Rotation: (b) Neutral to Maximum External Rotation: Total Included Angle: (b)	
of Motio	Internal/External Rotational Laxity Internal/External Rotational Laxity Internal/External Rotational Laxity Proximal/Distal Translational Laxity Internal/External Rotational Laxity		ning	Proximal/Distal Translational Laxity	Neutral to Maximum Proximal: (b)(4) Neutral to Maximum Distal: Total Proximal/Distal Laxity:
_			Pisto	Anterior/Posterior Translation Laxity	Neutral to Maximum Anterior: Neutral to Maximum Posterior: Total Anterior/Posterior Laxity:
				Neutral to Maximum Medial: Neutral to Maximum Lateral: Total Medial/Lateral Laxity:	
				All Humeral and Ulnar size combinations: (b)(4)	



Surgical Instrumentation Unique to the Device

Following are the instruments unique to the Zimmer Nexel Total Elbow. These are Class II accessories to the device. Applicable <u>material standards</u> are listed in a table at the end of this section.

	The Zimmer Nexel Fotal Elbow Class II Accessories to the Implant (Table Continued on Next Page)						
Catalog Number	Description	Material	Tissue Contact?	Image (not to scale)	Intended Use		
00-8401- 002-00	T-Handle	17-4 SS 303 SS 420 SS 316 SS	yes	Contraction of the second seco	Quick-connect handle for use with humeral and ulnar rasps.		
00-8401- 008-00	Reciprocating Saw Adapter	303 SS 420 SS 316 SS	yes	Contraction of the second seco	Adapter to connect humeral and ulnar rasps to a compatible large bonc reciprocating saw via the single use connector.		
00-8401- 008-01	Single Use Connector	455 SS	yes, skin only		Connects the reciprocating saw adaptor to a compatible large-bone reciprocating saw.		
00-8401- 009-00	Slide Hammer	17-4 SS 304 SS 17-4 SS	yes		Mates with humeral and ulnar stem extractors to deliver impact blows when extractor component is attached to the implant.		
00-8401- 018-00	Ulnar Bearing Tamp	17-4 SS	yes		To compress the articulating assembly (ulnar bearings and axle pin) into the humeral implant. (An alternate instrument to the articulation inserter.)		
00-8401- 012-04	Trephine Stabilizer, Size 4	17-4 SS	yes		Guides the finishing cut of the trephine through the		
00-8401- 012-05	Trephine Stabilizer, Size 5	17-4 SS	yes	↓ ST	anterior aspect of the humerus.		



	The Zimmer Nexel Total Elbow – Class II Accessories to the Implant (Table Continued from Previous Page, Continued on Next Page)					
Catalog Number	Description	Material	Tissue Contact?	Image (not to scale)	Intended Use	
00-8401- 014-07	Ulnar Provisional, Size 4, 75 mm, Left	17-4 SS	yes			
00-8401- 014-11	Ulnar Provisional, Size 4, 115 mm, Left	17-4 SS	yes			
00-8401- 015-07	Ulnar Provisional, Size 5, 75 mm, Left	17-4 SS	yes			
00-8401- 015-11	Ulnar Provisional, Size 5, 115 mm, Left	17-4 SS	yes			
00-8401- 016-09	Ulnar Provisional, Size 6, 90 mm, Left	17-4 SS	yes	(e		
00-8401- 016-11	Ulnar Provisional, Size 6, 115 mm, Left	17-4 SS	yes		To assess the bone preparation of the ulnar medullary canal, and to allow an intraoperative range of motion.	
00-8401- 024-07	Ulnar Provisional, Size 4, 75 mm, Right	17-4 SS	yes			
00-8401- 024-11	Ulnar Provisional, Size 4, 115 mm, Right	17-4 SS	yes			
00-8401- 025-07	Ulnar Provisional, Size 5, 75 mm, Right	17-4 SS	yes			
00-8401- 025-11	Ulnar Provisional, Size 5, 115 mm, Right	17-4 SS	yes			
00-8401- 026-09	Ulnar Provisional, Size 6, 90 mm, Right	17-4 SS	yes			
00-8401- 026-11	Ulnar Provisional, Size 6, 115 mm, Right	17-4 SS	yes			
00-8401- 019-00	Articulation Inserter	17-4 SS	yes		To compress the articulating assembly (ulnar bearings and axle pin) and ulnar implant into the humeral implant. Also, to compress the humeral bearing into the humeral implant. (This is an alternate instrument to the ulnar bearing tamp and the humeral bearing driver.)	
00-8401- 028-00	Ulnar Stem Inserter	Radel	yes		To transmit force to insert the ulnar implant without damaging the articular surface.	



77

Ĩ	The Zimmer Nexel Total Élbow – Class II Accessories to the Implant (Table Continued from Previous Page, Continued on Next Page)						
Catalog Number	Description	Material	Tissue Contact?	Imáge,(not to/scale)	Intended Use		
00-8401- 029-00	Implant Extractor Hook	17-4 SS	yes		To transmit force from the slide hammer to an implanted ulnar implant or humeral implant to facilitate their removal. The extractor hook is attached to the slide hammer by a threaded connection.		
00-8401- 033-01	Ulnar Rasp. Pilot, Left	17-4 SS	yes		Removes bone from		
00-8401- 033-02	Ulnar Rasp, Pilot, Right	17-4 SS	yes		the ulnar medullary canal and olecranon and shapes the medullary canal to		
00-8401- 034-01	Ulnar Rasp, Size 4/5, Left	17-4 SS	yes	~	create a precise geometry for a size- matched ulnar implant.		
00-8401- 034-02	Ulnar Rasp, Size 4/5, Right	17-4 SS	yes		Can be mated to the reciprocating saw adaptor or to the T-		
00-8401- 036-01	Ulnar Rasp, Size 6, Left	17-4 SS	yes		handle.		
00-8401- 036-02	Ulnar Rasp, Size 6, Right	17-4 SS	yes				
00-8401- 039-00	Ulnar Clearance Template	17-4 SS	yes		To check for adequate bone removal on the medial and lateral sides of the ulna for clearance for the articulation insertion instruments.		



	The Zimmer Nexel To (Table Continue			Accessories forfile É ontinued on Next Page}	mplant	
Catalog- Number	Description	Material	Tissue Contact?	Image (not to scale)	Intended Use	
00-8401- 044-10	Humeral Provisional, Size 4, 100 mm	17-4 SS	yes			
00-8401- 044-15	Humeral Provisional, Size 4, 150 mm	17-4 SS	yes			
00-8401- 045-10	Humeral Provisional, Size 5, 100 mm	17-4 SS	yes			
00-8401- 045-15	Humeral Provisional, Size 5, 150 mm	17-4 SS	yes		Dimensionally	
00-8401- 046-10	Humeral Provisional, Size 6, 100 mm	17-4 SS	yes		represents the humeral implant to help assess the preparation of the	
00-8401- 046-15	Humeral Provisional, Size 6, 150 mm	17-4 SS	yes		humeral medullary canal, and to intraoperatively check range of motion.	
00-8401- 055-15	Humeral Provisional, Long Flange, Size 5, 150mm	17-4 SS	yes			
00-8401- 055-20	Humeral Provisional, Long Flange, Size 5, 200mm	17-4 SS	yes			
00-8401- 056-15	Humeral Provisional, Long Flange, Size 6, 150mm	17-4 SS	yes			
00-8401- 056-20	Humeral Provisional, Long Flange, Size 6, 200mm	17-4 SS	yes			
00-8401- 058-04	Humeral Stem Inserter, Size 4	17-4 SS	yes		To insert the humeral implant into the humeral canal without	
00-8401- 058-05	Humeral Stem Inserter, Size 5/6	17-4 SS	yes	6.	damaging the humeral bearing.	
00-8401- 059-04	Humeral Extractor Plate, Size 4	17-4 SS	yes		To transmit force from slide hammer to humeral implant in order to facilitate	
00-8401- 059-05	Humeral Extractor Plate, Size 5/6	17-4 SS	yes		removal. Extractor plate is attached to humeral implant by the humeral extractor screws.	



	The Zimmer Nexel Fotal Elbow Class H Accessories to the Implant (Yable Continued from Previous Page, Continued on Next Page)						
Catalog Number	Description	Material	Tissue Contact?	Image (not to scale)	Intended Use		
00-8401- 060-00	Humeral Awl Reamer	17-4 SS	yes	×°	To open the humeral medullary canal. Also, provides a means of ensuring the trochlear resection is adequate to accept subsequent humeral rasps.		
00-8401- 061-00	I/E Alignment Rod	17-4 SS	yes		To assess and adjust, if necessary, the internal/external rotation of the humeral and ulnar rasps.		
00-8401- 064-00	Humeral Rasp. Pilot	17-4 SS	yes		Removes bone from the		
00-8401- 064-10	Humeral Rasp, Size 4, 100 mm	17-4 SS	yes		humeral medullary canal to create a precise geometry		
00-8401- 064-15	Humeral Rasp, Size 4, 150 mm	17-4 SS	yes		for a size-matched humeral implant. Humeral rasps also used to orient and securely		
00-8401- 065-10	Humeral Rasp, Size 5, 100 mm	17-4 SS	yes		support the humeral cut guide, and to guide the trephine cut.		
00-8401- 065-15	Humeral Rasp, Size 5, 150 mm	17-4 SS	yes		acpinic cut.		
00-8401- 066-10	Humeral Rasp. Size 6, 100 mm	17-4 SS	yes				
00-8401- 066-15	Humeral Rasp. Size 6, 150 mm	17-4 SS	yes				



	The Zimmer Nexel Total Elbow – Class II Accessories to the Implant (Table Continued from Previous Page, Continued on Next Page)						
Catalog Number	Description	Material	Tissue Contact?	Image (not to scale)	Intended Use		
00-8401- 070-45	Flexible Solid Reamer, 4.5 mm	13-8 SS	yes	ļ f	To remove bone from the ulnar medullary canal. Connects to a		
00-8401- 070-05	Flexible Solid Reamer, 5 mm	13-8 SS	yes		surgical drill or reamer driver.		
00-8401- 072-01	Flexible Cannulated Reamer, 5.5 mm	13-8 SS	yes				
00-8401- 072-02	Flexible Cannulated Reamer, 6 mm	13-8 SS	yes				
00-8401- 072-03	Flexible Cannulated Reamer, 6.5 mm	13-8 SS	yes] 			
00-8401- 072-04	Flexible Cannulated Reamer, 7 mm	13-8 SS	yes	ļ 📕	Mates over the guide wire. Used		
00-8401- 072-05	Flexible Cannulated Reamer, 8 mm	13-8 SS	yes		to remove bone from the ulnar and humeral medullary canals. Connects to a surgical drill or		
00-8401- 072-06	Flexible Cannulated Reamer, 9 mm	13-8 SS	yes	/	reamer driver.		
00-8401- 072-07	Flexible Cannulated Reamer, 10 mm	13-8 SS	yes				
00-8401- 072-08	Flexible Cannulated Reamer, 11 mm	13-8 SS	yes				
00-8401- 072-09	Flexible Cannulated Reamer, 12 mm	13-8 SS	yes				
00-8401- 074-00	Trephine, Size 4	17-4 SS 13-8 SS	yes		Removes bone from the distal humerus to accept a given size-		
00-8401- 075-00	Trephine, Size 5/6	17-4 SS 13-8 SS	yes		matched humeral implant. Connects to a surgical drill.		
00-8401- 078-00	Humeral Bearing Driver	17-4 SS	yes		To deliver a compressive force sufficient to fully seat the humeral bearing into the humeral implant when used with the humeral bearing driver pin. Also used to install a humeral bearing component into the humeral implant during a revision surgery. (This is an alternate instrument to the articulation inserter.)		



·	The Zimmer Nexel Total Elbów – Class II Accessories to the Implant (Table Continued from-Previous Page, Continued on Next Page)						
Catalog. Number	Description	Material	Tissue Contact?	Image (not to scale)	Intended Use		
00-8401- 079-00	Humeral Bearing Driver Pin	440A SS	yes		Subcomponent of the humeral bearing driver. Can also be used to aid in internal/external alignment of the ulnar provisionals.		
00-8401- 081-00	Ulnar Bearing Assembly Tool	17-4 SS 455 SS 316 SS	yes		To assemble the axle pin to ulnar bearings-B; and, to assemble the ulnar bearings-B onto the axle pin through the proximal head of the ulnar Implant. Performed by compressing the handle of the tool manually. To be used in situ.		
00-8401- 082-00	Humeral Bearing Placement Tool	420 SS 304 SS	yes	-	To place humeral bearing-A into the base of the humeral implant yoke.		
00-8401- 084-00	Humeral Screw Holder, Size 4	17-4 SS	yes		To securely hold the humeral screws and align the screws to a size-matched		
00-8401- 085-00	Humeral Screw Holder, Size 5/6	17-4 SS	yes		humeral implant.		
00-8401- 088-00	Humeral Screw Installation Driver Bit	440C SS LDPE	yes		Connects to the humeral screw installation driver to install the humeral screws into the humeral implant, and prevents the user from over-torquing the humeral screws. An alternative to the torque-limiting humeral screw driver, below. NOTE: <i>Cannot</i> mate with the humeral screw removal driver.		
00-8401- 089-00	Humeral Screw Removal Driver	455 SS 303 SS Silicone	yes		To remove the humeral screws from, the humeral implant only when mated with the humeral screw removal driver bit. NOTE: Will <i>not</i> allow installation torque to be applied.		
00-8401- 091-00	Humeral Bearing Extractor	17-4 SS	yes	Q	To remove the humeral bearing from the humeral implant.		
00-8401- 092-00	Articulation Extractor	17-4 SS	yes		To remove the articulation assembly (ulnar bearings and axle pin) and ulnar implant from the humeral implant.		



The Zimmer Nexel Total Elbow – Class II Accessories to the Implant (Table Continued from Previous Page, End of Table)					
Catalog Number	Description	Material	Tissue Contact?	lmage (not to scale)	Intended Use
00-8401- 093-00	Humeral Extractor Screw	455 SS	yes		Attaches the humeral extractor plate to the humeral implant to facilitate humeral component extraction.
00-8401- 097-00	Humeral Screw Removal Driver Bit	440C SS 316 SS	yes		Used with the humeral screw removal driver. NOTE: Cannot mate with the humeral screw installation driver.
00-8401- 004-00	Humeral Cut Guide, Size 4	17-4 SS	yes		Connects to humeral rasps. Used to facilitate the creation of parallel cuts
00-8401- 005-00	Humeral Cut Guide, Size 5/6	17-4 SS	yes	U	equidistant from the longitudinal axis of the rasp when used with a surgical blade.
00-8401- 080-00	Torque-Limiting Humeral Screw Driver	440C SST Ultem	yes		To install and remove the humeral screws into and from the humeral implant. Prevents the user from over- torquing the humeral screws upon insertion. This instrument is intended to be used for a single surgery only.



83

Surgical Instrumentation Not Unique to the Device

Following are the instruments not unique to the Zimmer Nexel Total Elbow. These are Class I instruments currently able to be used with other surgical devices cleared for market in the United States.

The Zimmer Nexel Total Elbow – Class I Instruments		
Catalog Number	Description	
00-4812-035-00	Small Hex Screwdriver	
47-2255-008-00	Sterile Ball Tip Guide Wire 2.4 x 70	
31-8106-168-00	Uinar Awi Reamer	
00-4811-035-01	Humeral Screw Installation Driver	
00-5049-053-00	Quick Use Curette	



84

Instrument Material Standards

The Zimmer Nexel Total Elbow: Instrument Material Standards					
Material	Applicable Standards per Zimmer Engineering Specifications				
6061-T6 Aluminum	ASTM B209, ASTM B211, ASTM B221				
13-8 SS	ASTM A484, ASTM A555, ASTM A564				
17-4 SS	ASTM A484, ASTM F899				
18-8 SS	ASTM A484, ASTM A555, ASTM A580				
303 SS	ASTM F899				
304 SS	ASTM A276, ASTM A484				
316 SS	ASTM F899				
316L SS	ASTM F138, ASTM F139, ASTM 1350, ISO 5832-1				
420 SS	ASTM F899				
440A SS	ASTM F899, ASTM A555				
440C SS	ASTM F899				
455 SS	ASTM F899, ISO 7153-1				
Phenolic	No Standard Applicable				
L 605 (Co-Cr-W-Ni alloy)	ASTM F90				
LDPE	ASTM D1248				
Radel	No Standard Applicable				
Silicone	No Standard Applicable				
Ultem	ASTM D5205				



85

Methods, Facilities and Controls

Method of Manufacturing





86

Packaging

Humeral Stem



The above listed materials are identical to those currently in use for Zimmer sterile implant packages with the only exception being that the configuration of the interior components correspond to the particular product being packaged.

Ulnar Stem



The above listed materials are identical to those currently in use for Zimmer sterile implant packages with the only exception being that the configuration of the interior components correspond to the particular product being packaged. Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017



Traditional 510(k) Premarket Notification

Articulation Kit



The above listed materials are identical to those currently in use for Zimmer sterile implant packages with the only exception being that the configuration of the interior components correspond to the particular product being packaged.

Humeral Screw Kit



The above listed materials are identical to those currently in use for Zimmer sterile implant packages. Only the configuration of the interior components changes to correspond to the particular product being packaged.

Pyrogenicity

This device is not labeled as nonpyrogenic. Per USP XXIII (161), requirements for specified endotoxin levels do not apply to orthopaedic implants.

Latex

There is no natural latex rubber in this product or its packaging.



Predicate Device(s)

The predicate and design basis for the Zimmer Nexel Total Elbow (the proposed device) is the Coonrad/Morrey Total Elbow, manufactured by Zimmer. Clearance letters (K001989 and K053189 cleared 25 July 2000 and 9 December 2005, respectively) for the predicate device are included at the end of this section.

Similarities and differences between the proposed implants and the predicate implants are tabulated below.

Substantial Equivalence Comparison

Indications for Use

	Proposed Device: The Zimmer Nexel Total Elbow	Predicate Device: The Coonrad/Morrey Totaf Elbow (K001918, K053189)
Intended Use	To be used in conjunction with bone cement to replace the articulating surfaces of the humeroulnar joint.	To be used in conjunction with bone cement to replace the articulating surfaces of the humeroulnar joint.
Indications for Use	 Elbow joint destruction which significantly compromises the activities of daily living Post-traumatic lesions or bone loss contributing to elbow instability Ankylosed joints, especially in cases of bilateral ankylosis from causes other than active sepsis Advanced rheumatoid, post-traumatic, or degenerative arthritis with incapacitating pain Instability or loss of motion when the degree of joint or soft tissue damage precludes reliable osteosynthesis Acute comminuted articular fracture of the elbow joint surfaces that precludes less radical procedures, including 13-C3 fractures of the distal humerus Revision arthroplasty CAUTION: This device is intended for cemented use only. 	 Elbow joint destruction which significantly compromises the activities of daily living Post-traumatic lesions or bone loss contributing to elbow instability Ankylosed joints, especially in cases of bilateral ankylosis from causes other than sepsis Advanced rheumatoid or degenerative arthritis with incapacitating pain Instability or loss of motion when the degree of joint damage precludes less radical procedures Revision arthroplasty CAUTION: This device is intended for cemented use only.
Product Code	ЛС	JDC




Technology

ł

	Feature	Proposed Device: The Zimmer Nexel Total Elbow	Predicate Device: The Coonrad/Morrey Total Elbow (K001918, K053189)
	Material	Wrought Tivanium Alloy (Ti-6Al-4V) with Titanium Plasma Spray	Wrought Tivanium Alloy (Ti-6Al-4V) with a porous coating of Titanium Beads
	Fixation Method	Intended for cemented use only	Intended for cemented use only
- General	Comparable Sizes	Size 4, Standard Flange, 100mm $\rightarrow \rightarrow \rightarrow$ Size 4, Standard Flange, 150mm $\rightarrow \rightarrow \rightarrow$ Size 5, Standard Flange, 100mm $\rightarrow \rightarrow \rightarrow$ Size 5, Standard Flange, 150mm $\rightarrow \rightarrow \rightarrow$ Size 5, Long Flange, 150mm $\rightarrow \rightarrow \rightarrow \rightarrow$ Size 6, Standard Flange, 100mm $\rightarrow \rightarrow \rightarrow$ Size 6, Standard Flange, 150mm $\rightarrow \rightarrow \rightarrow$ Size 6, Long Flange, 150mm $\rightarrow \rightarrow \rightarrow \rightarrow$ Size 6, Long Flange, 200mm $\rightarrow \rightarrow \rightarrow \rightarrow \rightarrow \rightarrow$	 → Size Extra Small, Standard Flange, 4in → Size Extra Small, Standard Flange, 6in → Size Small, Standard Flange, 4in → Size Small, Standard Flange, 6in → Size Small, Long Flange, 6in → Size Regular, Standard Flange, 4in → Size Regular, Long Flange, 6in
mponent	Sterility	Provided Sterile via Gamma Irradiation	Provided Sterile via Gamma Irradiation
Humeral Component – General	Design Fcatures	 Humeral stem designed to be implanted with bone cement into the patient's humeral medullary canal Anterior flange designed to accept a bone graft and limit torsional and posterior migration Humeral yoke with rounded corners to avoid the creation of stress risers within the medial and lateral humeral supracondylar columns Plasma spray region to enhance fixation to bone cement within the medullary canal, and to improve fatigue strength Features to accept Humeral Screws, Axle Pin and polymer Bearings A and B 	 Humeral stem designed to be implanted with bone cement into the patient's humeral medullary canal Anterior flange designed to accept a bone graft and limit torsional and posterior migration Humeral yoke with squared corners Titanium Bead coating region to enhance fixation to bone cement within the medullary canal Features to accept Axle Pin and polymer Bearings

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118






Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118 116



	Feature	Proposed Device: The Zimmer Nexel Total Elbow	Predicate Device: The Coonrad/Morrey Total Elbow (K004918, K053189)
	Material	Wrought Tivanium Alloy (Ti-6Al-4V) with Titanium Plasma Spray	Wrought Tivanium Alloy (Ti-6Al-4V) with Titanium Plasma Spray
	Fixation Method	Intended for cemented use only	Intended for cemented use only
Ulnar Component – General	Comparable Sizes	Left, Size 4, 115mm $\rightarrow \rightarrow \rightarrow \rightarrow \rightarrow \rightarrow \rightarrow \rightarrow \rightarrow$ Left, Size 5, 75mm $\rightarrow \rightarrow $	 → Left, Size Extra Small, Standard → Left, Size Extra Small, Long → Left, Size Small, Standard → Left, Size Regular, Standard → Left, Size Regular, Long → Right, Size Extra Small, Standard → Right, Size Extra Small, Long → Right, Size Small, Standard → Right, Size Small, Long → Right, Size Regular, Standard
IL Col	Sterility	Provided Sterile via Gamma Irradiation	Provided Sterile via Gamma Irradiation
Ulna	Design Features	 Left and right side-specific Ulnar stem designed to be implanted with bone cement into the patient's ulnar medullary canal Ulnar eye that is both highly polished and nitrogen-enriched to limit wear of the apposing polymer bearings Plasma spray region to enhance fixation to bone cement within the medullary canal 	 Left and right side-specific Ulnar stem designed to be implanted with bone cement into the patient's ulnar medullary canal Ulnar eye that is highly polished to limit wear of the apposing polymer bearings Plasma spray region to enhance fixation to bone cement within the medullary canal






- -



.

	Feature	Proposed Device: The Zimmer Nexel Lotal Fibow	Predicate Device: The Coonrad/Morrey Total Elbow (K001918, K053189)	
	Material	Vitamin E HXPE (Vitamin E Highly Cross-Linked Polyethylene)	UHMWPE (Ultra-High Molecular-Weight Polyethylene)	
-	Fixation Method	Press-Fit to Humeral Component	Sliding-Fit to Humeral Component Press-Fit to Ulnar Component	
inera	Sterility	Provided Sterile via Ethylene Oxide	Provided Sterile via Gamma Irradiation	
Polyethylene Components - General	Design Features	 Bearings A and B: Designed to broadly distribute joint reaction forces Designed to articulate against the Eye of the Ulnar Component Features to allow for press-fit mate to Humeral Component Bearings B only: Features to allow for press-fit mate to Axle Pin Designed to articulate against the inner-diameter and the medial and lateral surfaces of the Eye of the Ulnar Component Bearing A only: Designed to articulate against the outer-diameter of the Eye of the Ulnar Component 	 Humeral Bushings: Designed to articulate with the medial and lateral sides of the Eye of the Ulnar Component Features to allow for sliding-fit mate to the Humeral Component and sliding-fit mate to the Axle Pin (also known as Hinge Pin, also known as Snap Pin) Ulnar Bushing: Features to allow for a press-fit mate with the Eye of the Ulnar Component Designed to articulate against the Humeral Bushings 	







	Féature	Proposed Device: The Zimmer Nexet Total Elbow	Predicate Device; The Coonrad/Morrey Total Elbow (K001918, K053189)	
	Material	Zimaloy (Cobalt-Chromium-Molybdenum Alloy)	Zimaloy (Cobalt-Chromium-Molybdenum Alloy)	
Axle Pin	Fixation Method	 Press-Fit to Bearings B Compression fit between the taper of the Humeral Screws and the yoke of the Humeral Component 	 Sliding-Fit to Bearings Snap-Fit between inner and outer Axle Pin components 	
	Sizes	Size 4 (b)(4) Long Diameter • Solid Material Throughout Size 5/6 (b) Long (4) Diameter • Solid Material Throughout	Size Extra Small (b)(4) ' Long ' Outer Diameter ' Wall Thickness Size Small/Regular (b)(4) Long Outer Diameter Wall Thickness	
	Sterility	Provided Sterile via Ethylene Oxide	Provided Sterile via Gamma Irradiation	
	Design Features	 Designed to remain in a fixed position (not translate nor rotate) throughout the entire humeroulnar joint range of motion Machined complete as a single, solid component 	 Designed to freely rotate and slide Machined as 2 (an inner and an outer) hollow components that snap together intraoperatively 	
5	Material	Zimaloy (Cobalt-Chromium-Molybdenum Alloy)		
	Fixation Method	Threads engage with the Spiralock® threads of the Humeral Component		
	Sizes	Single size mates with all humeral component sizes and both axle pin sizes		
crews	Sterility	Provided Sterile via Gamma Irradiation		
Humeral Screws	Design Features	 Designed to remain in a fixed position throughout the entire humeroulnar joint range of motion (not directly loaded by any compressive joint reaction force during flexion or extension of the humeroulnar joint) Headless design intended to ensure the constant application of a clamp-load to the axle-pin and to avoid soft tissue impingement or irritation associated with a protruding screw head Non-threaded taper to (along with the Humeral Component) hold Axle Pin in compression Pilot feature to avoid cross-threading 	Not Applicable: No comparable part in the predicate system	



ł

l

		emi-Constraint Tharacteristics	Proposed Device: The Zimmer Nexel Total Elbow	Predicate Device: The Coonrad/Morrey Total Elbow (K001918, K053189)
mponent)	Varus/Valgus Angular Laxity		Neutral to Maximum Varus:Neutral to Maximum Varus:(b)((b)(Neutral to Maximum Valgus:(b)((c))((c))(Total Included Angle:(c))((b)((c))(
ained Range of Motion movement with respect to Humeral Component)		Internal/External Rotational Laxity	Neutral to Maximum Internal Rotation: (b) Neutral to Maximum External Rotation: (b) Total Included Angle: (b)	Neutral to Maximum Internal Rotation: (b)(Neutral to Maximum External Rotation: (b)(Total Included Angle: (b)
Motion	Pistoning	Proximal/Distal Translational Laxity	Neutral to Maximum Proximal: (b)(4) Neutral to Maximum Distal: Total Proximal/Distal Laxity:	Neutral to Maximum Proximal: (b)(4) Neutral to Maximum Distal: Total Proximal/Distal Laxity:
ange of] ent with	Pisto	Anterior/Posterior Translation Laxity	Neutral to Maximum Anterior: Neutral to Maximum Posterior: Total Anterior/Posterior Laxity:	Neutral to Maximum Anterior: Neutral to Maximum Posterior: Total Anterior/Posterior Laxity:
	-	Medial/Lateral Franslation Laxity	Neutral to Maximum Medial: Neutral to Maximum Lateral: Total Medial/Lateral Laxity:	Neutral to Maximum Medial: Neutral to Maximum Lateral: Total Medial/Lateral Laxity:
Semi-Constrained Range of Motion (defined in terms of Ulnar Component movement with respec			All Humeral and Ulnar combinations: (b)(4)	Humeral size XS with Ulnar size XS: (b)(4) Humeral size XS with Ulnar size S: (b)(4) Humeral size XS with Ulnar size REG: (b)(4) Humeral size S with Ulnar size XS: (b)(4) Humeral size S with Ulnar size S: (b)(4) Humeral size REG with Ulnar size XS: (b)(4) Humeral size REG with Ulnar size S: (b)(4) Humeral size REG with Ulnar size S: (b)(4) Humeral size REG with Ulnar size REG: (b)(4) Humeral size REG with Ulnar size REG: (b)(4)



Performance Specifications

	Test Results		
Test	Proposed Device: The Zimmer Nexel Total Elbow	Predicate Device: The Coonrad/Morrey Total Elbow (K001918, K053189)	
Fatigue Testing: Humeral Stem (Plasma Spray Region)	Five Zimmer Nexel Total Elbow humeral stems achieved (b)(4) load cycles at (b)(without fracture.	The predicate device humeral component was fatigue tested at various load levels, and the fatigue strength at ((b)(4) cycles was determined to be (4)(4)	
Fatigue Testing: Ulnar Stemstems achieved (h)(4)oad cycles in fatigue tesUlnar Stem (Mid-Stem Region)the mid-stem region at a load of (b)(without fracture.the fatigue cycles in t		The predicate device ulnar component was fatigue tested at various load levels, and the fatigue strength $a_1(h)(4)$ load cycles in the mid-stem region was determined to be (b)(4).	
Fatigue Testing: Ulnar Stem (Plasma Spray Region)	Five Zimmer Nexel Total Elbow ulnar stems achieved (b)(4)load cycles in the plasma spray region at a load of (b)(1) without fracture.	The predicate device ulnar component was fatigue tested at various load levels, and the fatigue strength at (4)(4) do ad cycles in the plasma spray region was determined to be (h)(
Wcar Testing	The Zimmer Nexel Total Elbow mean gravimetric wear rate $(h)(4)$ at $(b)(4)$ load cycles) was (h) less than that of the predicate.	The predicate device mean gravimetric wear rate was (h)(4)	
Durability Testing	The Zimmer Nexel Total Elbow device achieved (h)(4) run-outs at an equivalent of (h)(4) weight-in-hand.	The run out load for the predicate device was found to be equivalent to (b)(4) weight-in-hand.	
Modular Connection Fatigue Testing	The median fatigue strength of the non- articulating, mechanically locked, modular Zimmer Nexel Total Elbow implant components was no less than (b)(2)	Not Applicable: The Zimmer Coonrad/Morrey Total Elbow does not have non-articulating, mechanically-locked, modular implant components comparable to those of the proposed device.	

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service C-00412

JUL 2 5 2000

Ms. Laura D. Williams Regulatory Affairs Associate Zimmer P.O. Box 708 Warsaw, Indiana 46581-0708



Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Re: K001989 Trade Name: Coonrad/Morrey Total Elbow Regulatory Class: II Product Code: JDC Dated: June 28, 2000 Received: June 29, 2000

Dear Ms. Williams:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general control provisions of the Act. The general control provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the <u>Code of Federal Regulations</u>, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 - Ms Laura D. Williams

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for <u>in vitro</u> diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its Internet address "http://www.fda.gov/cdrh/dsmamain.html".

Sincerely yours,

ponne R. bochner.

Celia M. Witten, Ph.D., M.D. Director Division of General, Restorative and Neurological Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

Exhibit **B**

Page 1 of 1

510(k) Number (if known ____ KOO 1989

Device Name:

Coourad/Morrey Total Elbow

Indications for Use:

Indications include: post-traumatic lesions or bone loss contributing to elbow instability; ankylosed joints, especially in cases of bilateral ankylosis from causes other than sepsis; advanced rheumatoid or degenerative arthritis with incapacitating pain; revision arthroplasty, and instability or loss of motion when the degree of joint damage precludes less radical procedures.

The candidate for total elbow arthroplasty should exhibit joint destruction which significantly compromises the activities of daily living. Patients with single joint involvement (generally those with traumatic or degenerative arthritis) or significant lower extremity disability which require walking aids are less amenable to treatment than patients with advanced and predominately upper extremity involvement. If possible, elbow replacement should be done after hip or knee surgery to avoid excessive stress to the prosthesis required by crutch walking during total hip or knee rehabilitation.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use // (Per 21 CFR 801.109)

OR

Over-The-Counter Use (Optional Format 1-2-96)

(Division Sign-Off) Division of General Restorative Devices 510(k) Number <u>K 001989</u>

RA06004K.510


Records processed under FOIA request #2016-4653; Released by CDRH on 07/06/202005 08/1



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

DEC 9 2005



Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Laura D. Williams, RAC Manager, Corporate Regulatory Affairs Zimmer, Inc. P.O. Box 708 Warsaw, Indiana 46581-0708

Re: K053189

Trade/Device Name: Coonrad/Morrey Total Elbow Regulation Number: 21 CFR 888.3150 Regulation Name: Elbow joint metal/polymer constrained cemented prosthesis Regulatory Class: II Product Code: JDC Dated: November 11, 2005 Received: November 15, 2005

Dear Ms. Williams:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

-

Page 2 - Laura D. Williams, RAC

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <u>http://www.fda.gov/cdrh/industry/support/index.html</u>.

Sincerely yours,

Mark N. Melkerson Acting Director Division of General, Restorative and Neurological Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

16053189

Indications for Use

510(k) Number (if known):

Device Name:

Coonrad/Morrey Total Elbow

Indications for Use:

Indications include: post-traumatic lesions or bone loss contributing to elbow instability; ankylosed joints, especially in cases of bilateral ankylosis from causes other than sepsis; advanced rheumatoid or degenerative arthritis with incapacitating pain; revision arthroplasty, and instability or loss of motion when the degree of joint damage precludes less radical procedures.

The candidate for total elbow arthroplasty should exhibit joint destruction which significantly compromises the activities of daily living. Patient with single joint involvement (generally those with traumatic or degenerative arthritis) or significant lower extremity disability which require walking aids are less amenable to treatment than patients with advanced and predominantly upper extremity involvement. If possible, elbow replacement should be done after hip or knee surgery to avoid excessive stress to the prosthesis required by crutch walking during total hip or knee rehabilitation.

Prescription Use <u>X</u> (Part 21 CFR 801 Subpart D) AND/OR

Over-The-Counter Use _____ (21 CFR 807 Subpart C)

(Please do not write below this line - Continue on another page if needed)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off) Division of General, Restorative, and Neurological Devices


Page 1 of 1

510(k) Number_160531 84



Labeling

<u>Draft labels</u> for each of the proposed components are presented below, along with the <u>draft package insert</u> and <u>draft surgical technique</u>.



LOT Marker 99999999 ED: 008 Zimmer* Nexel™ Total Elbow Ulnar Component Plasma Sprayed Size 4 75 mm Length Left Tivanium® Ti-6AI-4V Alloy Sterile Qty-1 CCC Sterile Sterile Qty-1

EDI: 00840001407


REF Stars: 00-8400-014-07



Item Number: H_ProductID Description: H_ProductionDescript Operation: H_Operation Date of Issue: H_ProductionPackAppr Sterlity: H_Sterlitytype Prod Pack Version: H_ProductionPackVers Prod Pack Ver Date: H_ProductionPackDesi OtyIPkg. H_QuantityPerPackage

EDI NUMBER ver ediVer format ver formVer

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118



EDI: 00840004410 Zimmer⁴ Nexel¹⁴ Total Elbow Humeral Component Plasma Sprayed Size 4 100 mm Length Tivanium® Ti-6AI-4V Alloy Sterile Qty-1 REF 100-8400-044-10



Item Number, H_ProductID Description: H_ProductionDescript Operation: H_Operation Date of Issue: H_ProductionPackAppr Sterility: H_Sterilitytype Prod Pack Version: H_ProductionPackVers Prod Pack Ver Date: H_ProductionPackDesi QityIPkg: H_QuantityPerPackage

EDI NUMBER ver ediVer format ver formVer

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

132



LOT and 99999999 EDI: 008 Zimmer* Nexel** Total Elbow Humeral Component Plasma Sprayed Size 5 150 mm Length Long Flange Tivanium® Ti-6AI-4V Alloy Sterile Qty-1 CE (2) Automatic Rest

EDI: 00840005515

REF 100-8400-055-15



ttem Number: H_ProductID Description: H_ProductionDescript Operation: H_Operation Date of Issue: H_ProductionPackAppr Sterlinty: H_Sterlintytype Prod Pack Version: H_ProductionPackVers Prod Pack Ver Date: H_ProductionPackDesi OtyVPkg: H_QuantityPerPackage

EDI NUMBER ver ediVer format ver formVer

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

REF Hatter 00-8400-095-00 LOT Mumber 99999999 REF Heretence 00-8400-095-00 EDI: 00840009500 Zimmer * Nexel** Total Elbow Zimmer[®] Nexel[™] Total Elbow Articulation Kit Size 5,6 Articulation Kit Size 5.6 1 Axle Pin, 1 Humeral Bearing A, 2 Ulnar Bearings B Product Vivacit-E* Highly Crosslinked Polyethylene/Zimaloy* Co-Cr-Mo Alloy Sterile Oty-1 1 Axle Pin, 1 Humeral Bearing A, 2 Ulnar Bearings B Photo STERILE IEO Manufacturer Ster lized Using Instructions for Use ner Single Use oniv WARSAW IN 46580, U.S.A. 2019-09

Item Number: H_ProductID Description: H_ProductionDescript Operation: H_Operation Date of Issue: H_ProductionPackAppr Sterility: H_Sterilitytype Prod. Pack Version: H_ProductionPackVers Prod. Pack Ver Date: H_ProductionPackDesi OtyVPkg: H_OuantityPerPackage

EDI ver ediVer format ver formVer



Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

Floouce	REF Reference 00-8400-090-00 Zimmer® Nexel™ Total Elbow Humeral Screw Kit
Photo	

LOT 399999999	EDI: 00840009000	REF 100-8400-090-00
Zimmer* Nexel** Total Elba	ow'	
Humeral Screw Kit		
2 Humeral Screws		
Zimaloy® Co-Cr-Mo Alloy	Stenle Oty-1	We state
(CO A		TAT THE BODE REMANDO
Subalative Attertion, Ser	Stinker Using	WARSAW, IN 46580, U.S.A.
0058 shiple our planutatin lot by		

050

07/16/2017

Item Number: V Description: V Operation: V Date of Issue: V Sterility, V Prod. Pack Version: V Prod. Pack Ver, Date: V QtyVkg, V Records processed under FQIA request #2016-4653; Released by CDRH on 07/16/2017



English

.

The CE mark is valid only if it is also printed on the product label

ZIMMER® NEXEL™ TOTAL ELBOW

Records processed under FOIA request #2016-4653: Released by CDRH on 07/16/2017 service using a product placed on the market by Zimmer, the operating surgeon should study carefully the following recommendations, warmings and instructions, so well as the aveilable product-specific information (e.g., product literature, written surgical technique). Zimmer is not liable for complications arising from the use of this device outside of its indicated uses, not aligned with the surgical technique, product essection or similar matters outside the control of Zimmer.

DESCRIPTION

This device is a total elever positives designed for use with bone tement [] is available in sizes 4, 5 and 5, in left and right configurations. The ultrar and humeral impliant components are manufactured from *Trivanium*? Ti 6A1-4V alkey. The ultrar impliant component has a purcus coating of Ti-6A1-4V plasma spray and higs an anterior flange to accommodate a bone graft. The linkage pin and screws are present a program to solve a compared to account to be a compared to the program of a compared to compared to a compared to a compared to a compared to a co

INDICATIONS

- Indications for use include:
- Elbow joint destruction which significantly compromises the activities of daily tiving Post-traumatic lesions or bone loss contributing to elbow instability
- Ankylosed joints, especially in cases of bilateral ankylons from causes other than active sepsis
 Advanced meumatoid, post-traumatic, or degenerative arthntis with incapacitating pain
 Instability or loss of motion when the degree of joint or solt lissue damage precludes reliable
- osteosynthesis

.

- Acute communication and cutor fracture of the elbow joint surfaces that precludes less radical
 procedures, including 13-C3 fractures of the distal humanus
 Revision archroplasty

~ .

CAUTION: This device is intended for cemented use only

Records processed under FolA request #2016-4653; Released by CDRH on 07/16/2017

- Use of the Zimmer Nexer Total Elbow is contraindicated in patients with Currently active, or instory of repeated, local infection at the surgical site Paralysis or dysfunctional neuropathy involving the elbow joint

 - Significant resideral hand dysfunction Significant resideral hand dysfunction Excessive scaming of life skin or soft tissue that could prevent adequate soft lissue coverage Daily activities that would subject the device to significant stress (Fe , heavy labor, torsional stress,
 - and/or competitive sports)
 - Relative contraindications include
 - Distant fuci of infection (e.g. genifournary, pulmonary, sim (chronic lesions or ulcerationa), or other sites). In cases of distant infection, the fact of infection should be treated prior to, during and after
 - Surgery
 Ancient prior sepsis
 - WARNINGS

•

×.

- This device is for single patient use only. Do not reuse. Reuse of a single use device that has come in contact with blood, bone, has us only one tody fluids may lead to patient or user myre. Possible miss associated with reuse of a single use device include, but are not limited to, mechanical failure and transmission of infectious agents.
- To properly implicit this device, the sugreal lechnique should be consulted and carefully followed.
 Cerrenting technique is extremely important. The medullary canal should be consulted in the strenge blood fat, and bone debris and then thoroughly dred. A cement delivery system should be used to carefully fundaments.
- A snug mechanical fit of the implant within the humeral and uthar medullary canals may help minimize ÷. loosening of the implant components The amount of bone removed from the una should be sufficient to permit full elbow motion on the
- ς. operating table when the appropriate provisional positres is sfully inserted. The surgeon should check for full range of elbow motion at appropriate times during surgery.

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017

- Proper handling of this implant is important. Contouring (bending) of the humeral or ultrar stems should be avoided. Any alteration may produce defects that could become focal points for excessive stresses leading to implant failure.
- To prevent late infection of total york implants many surgeons advise the use of antibiotic prophylaxis before and after dental manipulation, endoscopic examinations and other minor surgical procedures. ε.
- Patents with significant lower extremity disability who require waking acts are less amenable to total elow arthroplasty inan patients with advanced and predominately upper extremity involvement.

ADVERSE EFFECTS

- The following adverse effects have been reported.
- · Early or late component loosening Implant failure or fracture
- Infection
- Nerve injury
- Triceos avulsion
- Cortical perforations and fractures
- Adverse tasue response and inflammatory reaction to the implant material or wear debits, including • metal sensitivity
- Osteolysis and/or loosening due to wear debns.

STERILITY

STERULTY The Zimmer Nexe/ Total Elbow implants are provided stenie (the humeral and ulmar implant components, and the humeral screws are stenized by gamma irrediation, indicated by the <u>STERUE</u> <u>R</u> symbol on the labeling, and, the articulation kit – including the bearings and the axie pin – is stenized by athylane oxide, indicated by the <u>STERUE</u> <u>SO</u> symbol on the labeling) and remain stenie until the expiration date displayed on the package, as long as the package inlegity has not been indicated inspect each package prior to use and do not use the implant component if any seal or cavity is damaged or breached on if the expiration date has been exceeded. Once opened, the implant component must be und displayed to an exceeded. used, discarded or re-stenlized

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017

These startization instructions are consistent with ANSI/AAMI/ISO standards and guidelines. They should be used for startie items that were opened but unused Do not reuse implants.

- Do not reuse implants

.

- In the event of inadvertent loss of stanlity while preparing for surgery, all metal implant components may be re-stanlized only once for immediate use. This is subject to the exceptions listed below <u>DO NOT RE-STERILIZE</u>

 - Impaint components that have been contaminated with body fluids or biologic debns, or that have been previously implanted
 - Implant components on the accuracy implantation date that has been exceeded Implant components containing Weard-E polyethylene
- Impaint comparisits containing views to possible transmission or possible real of the provided of the provided of the provided of the table under forms to the recommended specifications for stars technization provided in the table under Recommended Starkzahor/Krestaritization Specifications before the baseline of the table under starge enough to contain the employed the provided of the table under starge enough to contain the employed provided the baseline of the provided starkzahor/Krestaritization Specifications before the provided starkzahor/Krestaritization starkzahor/Krestaritization starkzahor/Krestaritizations before the provided starkzahor/Krestaritization starkzahor/Krestaritization starkzahor/Krestaritization starkzahor/Krestaritization starkzahor/Krestaritization starkzahor/Krestaritization starkzahor/Krestaritization starkzahor/Krestaritization starkzahor/Krestaritization starkzahor/Kres
- implant component without stressing its seals or learing Implant components must be disassembled and sterilized separately to minimize potential bio-burden buildup in the dead space and expansion/contraction streases.
- Risse implant components, paying special attention to porous coatings, to remove lint or other non-biologic debris using USP (United States Pharmacope a) purfied water .

140

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017 Follow the stenizer manufacturer's instructions by partients and selection of storkization parameters. Drying littles vary according to load size and should be increased for larger loads.

For Single, All Metal Implants:

	Stearn Sten	lization		
Туре	Temperature	Exposure Time	Mmimum Dry Time	
Gravity Displacement	121°C (250°F)	30 minutes		
Gravity Displacement	132°C (270°F)	15 minutes	_	
UK Pre-vacuum/Pulsaling Vacuum ¹	134°C (273°F)	3 minutes	15 minutes	
Pre-vacuum/Puisating Vacuum	132°C (270°F)	4 minutes		

* This cycle is not for use in the United States

Please contact Zimmer at the following number if you have additional questions. In the USA, call 1-800-348-2759 For calls outside the USA, call the local international access code +1-574-267 6131

.

- STERILIZATION PRECAUTIONS Aggressive cleaning with delergents and brushes may damage special features of the implants, such as portula coatings. Also, defain detergents may be difficult to ninse off polymeritems: especially those made of silicone rubber
- Items made from intanium and bitanium alloys can form oxide layers from steam borier treatment chemicals or detergent residues. While these paces are traccompatible, they can obliterate identification etchings and stampings.

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017

Complications and/or failure of possibility mplants are more likely to occur in patients with unrealistic functional expectations, heavy pagents, physically active patients, end/or with patients who fail to follow through with the required rehabilitation program. Physical activity can result in loosening, wear, and/or fracture of the implant. The prospective implant patient must be counseled about the capabilities of the implant and the impact it will have on his or her lifestyle. The patient must be instructed about all it postoperative restrictions particularly those related to occupational and sports activities, and about the possibility that the implant or its components may wear out, fail or need to be replaced. The implant may not last the rest of the patient's life, or any particular length of time. Because prosthetic implants are not as strong, reliable, or durable as natural, healthy tissues/bones, all such devices may need to be replaced at some point.

Note. The prosthasis will not restore function to the level expected with a normal healthy joint and has functional irmitations. Excessive muscular activity such as pounding and or carrying excessive loads must be evolved or restricted. The patient must not 1 if more than one pound (~0 kg) during the first three postoperative months, and, thereafter, not more than five pounds (~2.25 kg) with the operated arm in addition, the nake and complications listed under ADVERSE EFFECTS must be explained to and discussed with the patient.

CAUTION: This device is intended for comented use only

. .

Draft Surgical Technique

Preoperative Considerations Indications/Contraindications Short Version – Reference Surgical Approach Humeral Preparation Ulnar Preparation Trial Reduction Bone Graft Preparation Ulnar Component Implantation Humeral Component Implantation Closure Postoperative Management Revision Procedure Scrub Nurse Instructions

Pre-Operative Considerations

For those inexperienced in the technique of elbow arthroplasty, training with a cadaver specimen(s) is recommended

The surgeon should:

- Be aware of how to articulate and disarticulate the new coupling mechanism
- Understand the technique to place a bone graft beneath the anterior flange of the humeral component
- Be careful of implant sizing and cement technique in patients with pre-existing shoulder pathology
- Understand if a long humeral stem is to be implanted, stem modification may be needed if significant anterior bowing is noted on the lateral pre-operative radiographs
- Surgeon should consider counter-sinking the **Humeral Component** up to 10 mm proximal to the axis of rotation to help facilitate post-op range of motion
- The amount of bone removed from the ulna (including the tip of the olecranon and coronoid process) should be sufficient to permit full elbow motion on the operating table when the appropriate Humeral and Ulnar Provisional prostheses are fully inserted. The surgeon should check for full range of elbow motion at appropriate times during the surgery
- Be aware of pinch-points, sharp edges and secure handling of components (i.e. transferring components from back table to operating table) while assembling implants and utilizing instruments
- Be aware that size "5/6" instruments indicate that this instrument can be used for a Size 5 or 6 implant; likewise "4/5" indicates applicability to a Size 4 or 5 implant

Indications / Contraindications

INDICATIONS

Indications for use include:

- Elbow joint destruction which significantly compromises the activities of daily living
- Post-traumatic lesions or bone loss contributing to elbow instability
- Ankylosed joints, especially in cases of bilateral ankylosis from causes other than active sepsis
- Advanced rheumatoid, post-traumatic, or degenerative arthritis with incapacitating pain
- Instability or loss of motion when the degree of joint or soft tissue damage precludes reliable osteosynthesis
- Acute comminuted articular fracture of the elbow joint surfaces that precludes less radical procedures, including 13-C3 fractures of the distal humerus
- Revision arthroplasty

CAUTION: This device is intended for cemented use only

CONTRAINDICATIONS

Use of the Zimmer Nexel Total Elbow is contraindicated in patients with:

- Currently active, or history of repeated, local infection at the surgical site
- Paralysis or dysfunctional neuropathy involving the elbow joint
- Significant ipsilateral hand dysfunction
- Excessive scarring of the skin or soft tissue that could prevent adequate soft tissue coverage
- Daily activities that would subject the device to significant stress (i.e., heavy labor, torsional stress, and/or competitive sports)

Relative contraindications include:

- Distant foci of infection (e.g. genitourinary, pulmonary, skin [chronic lesions or ulcerations], or other sites). In cases of distant infection, the foci of infection should be treated prior to, during and after surgery.
- Ancient prior sepsis

WARNINGS

- This device is for single patient use only. Do not reuse. Reuse of a single use device that has come in contact with blood, bone, tissue or other body fluids may lead to patient or user injury. Possible risks associated with reuse of a single use device include, but are not limited to, mechanical failure and transmission of infectious agents.
- To properly implant this device, the surgical technique should be consulted and carefully followed.
- Cementing technique is extremely important. The medullary canal should be copiously irrigated to remove blood, fat, and bone debris, and then thoroughly dried. A cement delivery system should be used to cement the humeral and ulnar implant components.
- A snug mechanical fit of the implant within the humeral and ulnar medullary canals may help minimize loosening of the implant components.
- The amount of bone removed from the ulna should be sufficient to permit full elbow motion on the operating table when the appropriate provisional prosthesis is fully inserted. The surgeon should check for full range of elbow motion at appropriate times during surgery.

.

.

.

Short Version – Reference

-

3

Include a page of mostly images (cliff notes version of technique right up front)

Chapter 1: Surgical Approach


1.0 Tourniquet

• Apply the tourniquet

1.1 Incision

- Position
 - The recommended position is supine, with a sandbag or rolled towel under the scapula and arm placed across the chest
- Initial incision
 - Make a straight incision over the elbow joint, approximately 15cm in length
 - Center the incision just lateral to the tip of the olecranon (Figure 1.1)

Note: a more mid-line positioned incision decreases the need for elevating extensive flap






1.2 Protect the Ulnar Nerve

- Identify the ulnar nerve
 - Identify the medial aspect of the triceps mechanism and isolate the ulnar nerve using ocular magnification
- Transpose the ulnar nerve
 - Mobilize the ulnar nerve well proximally and distally to the first motor branch and carefully translocate it anteriorly into the subcutaneous tissue (Figure 1.2)
 - o Utilize the bipolar cautery for hemostasis
 - o Incise or excise a small portion of the medial intramuscular septum

Note: The nerve must be carefully protected throughout the remainder of the procedure



The Bryan-Morrey Approach

The standard surgical approach employs Bryan/Morrey exposure and is recommended for new users of the New Zimmer Total Elbow system and users inexperienced with power rasping. The Bryan/Morrey approach employs a meticulous repair of the triceps that is detailed at the end of this surgical technique. Once experience is gained, other exposures (e.g. Triceps-On/Sparing and Triceps-Split) and/or power-rasping of the medullary canals can be employed at the surgeon's discretion.

- Release the triceps
 - Make an incision just medial to the crest of the ulna and elevate the ulnar periosteum along with the triceps and forearm fascia (Figure 1.3A)
 - Dissection proceeds laterally insuring that the triceps remains in continuity with the periosteum and fascia
 - Sharpey's fibers are carefully released and the anconeus elevated from its ulnar insertion




Figure 1.3A

- Visualize the joint
 - Dislocate the entire extensor mechanism (triceps, ulnar periosteum, and anconeus) as a single soft-tissue sleeve laterally
 - This allows complete exposure of the posterior aspect of the distal humerus, proximal ulna, and radial head (Figure 1.3B)
 - Release the anterior capsule and contracted soft tissue from the distal humerus to eliminate any flexion contracture. (Figure 1.3C)



Figure 1.3B



Figure 1.3C

- Disarticulate the ulnohumeral joint
 - Subcutaneously release the medial and lateral collateral ligaments from their humeral attachment
 - Flexion of the elbow to disarticulate the ulna from the humerus (Figure 1.3D)
 - If further visualization is required, external rotation of the forearm can facilitate exposure of contracted soft-tissue

Note: A complete release of the soft tissues around the medial side of the distal humerus is intended to protect the medial epicondyle (the weakest part of the distal humerus) from fracture during flexion and manipulation of the forearm. Elbows with severe arthritis, post-traumatic surgery, and/or asymptomatic soft-tissue contractures should have the entire capsule and soft tissues from the distal humerus circumferentially released subcutaneously.





Figure 1.3D

Chapter 2: Humeral Preparation

2.0 Resect the Trochlea

- Visualize the width and location of the initial resection
 - Hold the barrel of the central portion of the **Humeral Awl Reamer** against the trochlea and visually align the long axis of the instrument with the long axis of the humerus (Figure 2.0A)
- (Scrub Nurse will load appropriately sized saw blade)
- Resect approximately 20 mm of the middle portion of the trochlea using a saw or rongeur (Figure 2.0B)
 - The width of the resection should be sufficient to avoid impingement of the **Humeral Rasp**
 - Retain trochlear bone for anterior distal bone graft (Chapter 5)

Tip: The bone graft, which measures 3-5mm thick, can be cut by first making a center cut, followed by additional oblique medial & lateral cuts. One of these slices can be used as the anterior flange bone graft



2.1 Expose the Humeral Medullary Canal

- Use a bur or rongeur to penetrate a small portion of cortical bone from the proximal roof of the olecranon fossa (Figure 2.1)
- Proceed with control until cancellous bone can be seen



Figure 2.1

2.2 Ream the Humeral Medullary Canal

- Insert the shaft of the **Humeral Awl Reamer** into the humeral medullary canal and lightly twist in order to open the canal and verify the axis (Figure 2.2)
 - The barrel of the **Humeral Awl Reamer** should be centered and fit through the previous trochlear cut, otherwise remove more bone until it fits



Note: the barrel on the central portion of the **Humeral Awl Reamer** ensures clearance for the width of the largest **Humeral Rasp**

2.3 Humeral Canal Preparation

Iscrub Nurse section will describe how to attach Rasps)

Note: The posterior face of the **Rasp** should be kept parallel to the plane formed by the posterior cortices of the medial and lateral columns. If these landmarks are not available, the relatively flat posterior surface of the shaft of the distal humerus can also be used to approximate this plane (Figure 2.3A)

Note: The orientation of the **T-Handle** can be altered 90 degrees to suit surgeon preference





Figure 2.3A&B

- Advance the **Pilot Humeral Rasp** manually until it meets resistance, then gently impact it proximally into the canal until the solid etched line on the **Rasp** aligns with the anatomic axis of flexion (Figure 2.3B) (*see pre-op considerations*)
- Serial rasping continues until the **Rasp** fits snugly in the canal at a depth at which the flexion axis is in alignment with the solid etched line on the **Rasp**

- Confirm Humeral Rasp alignment by inserting the larger diameter end of the I/E
 Alignment Rod through the Humeral Rasp to verify it is in same plane as the posterior cortices or the flat posterior surface of the distal humerus if cortices not present (Figure 2.3C)
- Do not remove the final Humeral Rasp from the canal as it is used to direct the trephine cut below

Note: The **Rasps** are marked to indicate implant features (Figure 2.3C)

Tip: Flexible Cannulated Reamers are available to safely enlarge the canal if even the smallest Rasp is too large for the canal




Figure 2.3C

2.4 Preliminary Trephine Cut

- Choose the appropriately sized **Trephine**, either a Size 4 or 5/6 based on the final **Rasp** utilized
- Insert its pilot pin of the Trephine into the hole of the Humeral Rasp (Figure 2.4A)

Tip: the use of a high speed attachment is recommended

...




Figure 2.4A

Note: Size 4 Trephine has a slightly different diameter than the Size 5/6 Trephine to prevent mismatch from occurring between the Trephine and Humeral Rasp

- Start **Trephine** drill before its teeth are engaged in the posterior humeral bone, then gently advance until the pilot pin reaches its depth stop (Figure 2.4B)
 - This will score the posterior surface and serve as a guide for the final **Trephine** cut. Leave bone fragments in place



2.5 Make Vertical Cuts

- Attach the appropriate Humeral Cut Guide to the Humeral Rasp; choose Size 4 or 5/6 (Note: if depth of rasp insertion adjustments are desired, be sure to impact the T-Handle and not the Cut Guide)
- Stabilize the Humeral Cut Guide with the Humeral Bearing Driver Pin (Figure 2.5A&B)



- Using the slots of the **Humeral Cut Guide**, make required medial and lateral cuts; this determines the final width of the humeral preparation (Figure 2.5B)
- Remove the instruments

Note: for the final version of the surgical technique, this image will show the Humeral Cut Guide and Rasp and Saw Blade, with humerus. It will also show retractors to protect soft tissues and to avoid saw blade.




Figure 2.5B

2.6 Finish the Trephine Cut


- Insert the appropriately sized **Trephine Stabilizer** into the humeral medullary canal and align the center of the **Trephine Stabilizer** hole with the center of previous trephine cut; choose from Size 4 or 5/6 (Figure 2.6A)
- Drill until Trephine pilot pin reaches the depth stop
- Remove instruments and bone fragments (Figure 2.6B)





2.7 Assess depth and orientation of the Humeral Provisional

- Do not cement
- Place the appropriate **Humeral Provisional** into the prepared humeral canal (Figure 2.7A)
 - Ensure that the **Humeral Provisional** is fully seated (impaction should not be required), the flexion axis is in the proper position, and that it does not sit proud relative to the distal aspect of the lateral humeral condyle
- Trim any excess condylar bone distal to the provisional with rongeurs
- Remove instruments



-

Chapter 3: Ulnar Preparation

3.0 Ulnar Canal Preparation

- Remove the tip of the olecranon (Figure 3.0A) and use a high-speed bur to penetrate • cortical bone at the base of the coronoid to enter the intramedullary canal (Figure 3.0B)
- Create a "notch" along the guiding ridge of the ulna. This is a critical step to avoid placement of the ulnar component in a flexed position
- Place a finger over the exposed proximal shaft of the ulna to help prevent violation of the medullary canal in anticipation of the ulnar bow (Figure 3.0C)
- Carefully advance the Ulnar Awl Reamer axially into the medullary canal (Figure 3.0D)



Figure 3.0A





Figure 3.0B



Figure 3.0C

Note: for the final version of the surgical technique, this image will be updated to show in situ




3.1 Ulnar Canal Reaming (if necessary)

Note: the 4.5mm **Reamer** is prescribed for the size 4 implant; the 6.5 **Reamer** for the size 5, and; the 7.0 **Reamer** for the size 6 ulnar devices.


Reamers must be used in progressively larger cutting head sizes. Always begin with the 4.5mm solid reamer. DO NOT skip sizes, or attempt to begin with larger cutting head sizes.

• Advance the 4.5 mm Flexible Solid Reamer carefully until the etch line of the desired Ulnar Component length (75, 90, or 115mm) reaches the entrance to the ulnar canal (distal to the greater sigmoid notch) (Figure 3.1A)



- Continue to ream canal sequentially
 - Flexible Solid Reamers are size 4.5mm and 5mm, while Flexible Cannulated Reamers are sizes 5.5 7.0mm (Figure 3.1B)

Note: Flexible Cannulated Reamers require use with a **Sterile Ball Tip Guide Wire 2.4 x 70** to ensure proper tracking and control depth of insertion



3.2 Rasp the proximal ulnar canal

- Attach **T-Handle** to the **Pilot Ulnar Rasp** and properly orient the **Rasp** in both the coronal and sagittal planes; choose from left or right **Ulnar Rasp** configurations (Figure 3.2A)
- Confirm Ulnar Rasp alignment by inserting the narrow end of the I/E Alignment Rod through the Ulnar Rasp to verify it is parallel to the relatively flat surface of the posterior aspect of the Olecranon (Figure 3.2A)
- Advance the **Pilot Ulnar Rasp** manually until it meets resistance, then gently impact it distally into the canal until the center of the Ulnar eye is aligned to the projected center of the greater sigmoid notch
- Continue with sequential rasping until the proper size has been determined (Size 4/5 **Rasp** is the final **Rasp** used for a Size 4 or 5 implant; Size 6 **Rasp** is the final **Rasp** used for a Size 6 implant)
- Leave the final **Rasp** in canal




Figure 3.2A

Note: Features and markings on the medial, lateral and posterior faces of the **Rasps** correspond to characteristics of the **Ulnar Component** (Figure 3.2B)

Image show Ulnar Component along side ulnar rasp demonstrating the feature similarities

Figure 3.2B

3.3 Sigmoid Notch Preparation

- Place Ulnar Clearance Template through the Ulnar eye feature of the Rasp from the medial and lateral side in order to ensure adequate clearance for Ulnar Bearings and final assembly step (Figure 3.3)
- Remove any impinging bone to allow for proper seating of the Humeral Component



rigure 5.5

3.4 Assess the accuracy of the ulnar canal preparation

- Do not cement
- Place Ulnar Provisional into the ulnar canal (impaction should not be required) and confirm that the center of the circular head of the Ulnar Provisional is concentric with the projected center of the greater sigmoid notch (Figure 3.4A)
 - Humeral Bearing Driver Pin can be used to assess proper internal/external alignment (Figure 3.4B)
- Do not remove the Ulnar Provisional once proper seating is established

Image will show ulnar provisional in ulna, lateral view and the distal (transverse view) showing into of humeral bearing driver pin sliding into through hole

Figure 3.4A&B

Chapter 4: Trial Reduction

- Do not cement
- Re-insert the appropriate Humeral Provisional into the humerus
- Slide the Ulnar Provisional into the Humeral Provisional (Figure 4.0)
- Perform a trial range of motion
- Remove any bone impingements (possibly including portions of the radial head) with a bur or rongeur and perform any soft tissue release as needed

Note: **Provisionals** will not lock together and will not mate with implants. Causes for incomplete restoration of elbow extension include too distal orientation of humeral joint axis, inadequate depth of insertion of ulnar provisional, unresolved angular deformity, inadequate release of anterior soft-tissue contracture and posterior bone impingement



Figure 4.0

Chapter 5: Bone Graft Preparation

- A bone graft is fashioned from an excised trochlea or radial head
- Trim graft to size with an oscillating saw or rongeur
 - Approximately 2-3mm thick, 2 cm long and 1-1.5mm wide (Figure 5.0)

Bone graft			
E: E 0			



Chapter 6: Ulnar Component Implantation

- Clean the ulnar canal with copious irrigation to remove blood, fat and bone debris, and then thoroughly dry.
- Cut the cement tube to the length of the Ulnar Component (Figure 6.0A)

- Use a retrograde technique to fill canal with cement
 - Leave approximately 1 cm of the proximal canal free of cement to avoid excessive backflow

Tip: to help avoid porosity in the cement mantle, use finger to compress the cement into the canal (gauze will protect glove from cement.) Trim the cement tube flush with the nozzle, apply the nozzle to the canal and continue injecting cement until cement visibly backflows around nozzle, clear the excess cement




Figure 6.0A

- Insert the Ulnar Component so the center of Ulnar eye is aligned with the projected center of the greater sigmoid notch (Figure 6.0B)
 - o Carefully impact Ulnar Component with the Ulnar Stem Inserter to fully seat
 - Use the Quick Use Curette to thoroughly remove any extruded bone cement displaced by the stem insertion
- DO NOT install the Articulation Assembly (consisting of the Ulna Component, Axle-Pin and two Ulnar Bearings) until after the Ulnar Component has been placed properly in the canal and all bone cement has been removed from the exposed articulation area




Figure 6.0B

Chapter 7: Humeral Component Implantation

Note: The **Humeral Component** is cemented. The plasma spray yoke and anterior flange is for bone growth. No cement is utilized on these surfaces.

- (Scrub Nurse will install the Humeral Bearing into the Humeral Component)
- Prep the canal by using copious irrigation to remove blood, fat and bone debris, and then thoroughly dry
- Place cement restrictor in place
- Inject bone cement using a retrograde fill technique
 - Cut the cement tube to the length of the Humeral Component (Figure 7.1A)

• Apply cement to a depth 2 cm distal to the tip of the length of the **Humeral Component**; leave approximately 1 cm of the distal canal free of cement to avoid excessive backflow

Alternative Method: to help avoid porosity in the cement mantle, use finger to compress the cement into the canal (gauze will protect glove from cement.) Trim the cement tube flush with the nozzle, apply the nozzle to the canal and continue injecting cement until cement visibly backflows around nozzle, clear the excess cement




Figure 7.1A

- Insert the Humeral Component into the medullary canal
 - Wedge the bone graft between the flange of the **Humeral Component** and the anterior distal humeral cortex
 - Ensure the graft is engaged by the Humeral flange during insertion
- Carefully impact the Humeral Component securely in place with the Humeral Stem Inserter (Figure 7.1B)
- Clear any excess bone cement
 - Pay special attention to remove cement from the posterior hole on the humerus that is used to remove the humeral bearing during revision cases




Note: for the final version of the surgical technique, this image will show the humeral component with preloaded Humeral Bearing. It will also label the Humeral stem inserter.

Figure 7.1B

7.2 Attach the Bearings/Axle-Pin to Ulnar Component


- (Scrub Nurse will prepare Bearing Assembly at back-table)
- Perform installation in situ, after the Ulnar Component has been cemented in place

- Introduce the prepared Ulnar Bearings/Axle-Pin Assembly to the implanted Ulnar Component
 - Use caution to avoid contact between the Axle-Pin and the Ulnar eye to avoid scratching the eye's surface (Figure 7.2A)



• Squeeze the handles of the Ulnar Bearing Assembly Tool until the Axle-Pin has been fully seated into the opposing Ulnar Bearing and resistance is felt


Note: A visual cue on the handle will confirm when the stopping point for Size 4 and 5/6 has been met (Figure 7.2B)




7.3 Reduce the Elbow

- Allow cement to cure
- Apply the feet of the Articulation Inserter to the pocketed tabs of the Ulnar Bearings and partially reduce the joint by guiding the ends of the Axle-Pin into the Humeral yoke (Figure 7.3A)
- Line up the Ulnar Bearing tabs with the Humeral Component slots
- While maintaining the Articulation Inserter's feet on the Ulnar Bearings, insert the peg of the proximal leg of the Articulation Inserter into the hole at the base of the Humeral yoke (Figure 7.3B)
- Squeeze the handle of the Articulation Inserter to drive the Ulnar Bearings into the Humeral yoke, ensuring that the tabs of the Ulnar Bearings maintain alignment with the slots of the Humeral yoke throughout insertion (Figure 7.3C)
- Remove tools

Note: The Ulnar Bearing Tamp is an alternate tool available to assist with this step



- Avoid contact of the Ulnar eye with the Articulation Inserter
- Full insertion is achieved when the surfaces of the Ulnar Bearings appear to be flush with the distal Humeral yoke (Figure 7.3E)



Note - the following edits will be made for the final version of the surgical technique: show inset of the screw hole ghosted and the position of the Axle pin to demonstrate how the screws ride above axle.

7.4 Insert Humeral Screws

Do not tighten screws with an alternate driver, do not use power. Use the tools provided in the instrument set to insert screws. Proper application of torque to install the Humeral Screws is required for a successful implant.

- E (Scrub Nurse will prepare the Screws/Screw Holder and Driver Bit/Driver)
- Insert each screws from the Humeral Screw Holder and into the Humeral Component yoke (Figure 7.4A)
 - o Turn the Driver until the first screw is free of the instrument (7.4B)
 - Repeat process with the second screw
 - Tighten the screws in two sequences (screw 1 then screw 2, followed by screw 1 and screw 2 again) with the **Driver** until an abrupt loss of torque is felt
 - The **Driver Bit** is designed to yield when the appropriate torque is applied to the screw, utilizing a protective sleeve to keep both parts of the yielded bit contained
 - o The Driver Bits are designed and calibrated for single-use/single-surgery

Note: An optional Elbow Torque Driver may be used, which will produce an audible "click" once torque has been achieved; Elbow Torque Driver is designed for single-use/single-surgery




Figure 7.4B

7.5 Perform final range of motion

- Perform a final range of motion
- Remove any bone impingements (possibly including portions of the radial head) with a bur or rongeur and perform any soft tissue release as needed

Chapter 8: Closure

8.0 Close the Wound

- Repair the Triceps (Figures 8.0A-E)
 - o Slightly over correct the triceps mechanism to its anatomic position
 - o Place cruciate and transverse holes in the proximal ulna
 - Use a heavy #5 nonresorbable suture in a crisscross fashion
 - A transverse suture compresses the triceps against the olecranon
 - Tie sutures with the elbow in approximately 45 degrees of flexion
- Complete the closure in a routine fashion
 - Stabilize the ulnar nerve in the subcutaneous anterior pocket
 - Deflate the tourniquet and obtain hemostasis
 - Close the wound in layers
 - o Insert a drain, if desired
- Finish closure
 - \circ Apply a compressive dressing and split with the elbow in full extension





Figure 8.0B



Figure 8.0C



Chapter 9: Postoperative Management

9.1 Postoperative Management

- Apply a light dressing and elevate arm
- 24 hours
 - o Remove the drain, if used, at approximately 24 hours
- 36-48 hours
 - Remove the compressive dressing on the second day after surgery
 - Instruct the patient on activities of daily living
 - Typically, no formal physical therapy is required or indicated unless necessary for the shoulder or hand
 - Avoid strengthening exercises
 - Allow elbow flexion and extension as tolerated
- 1-3 Days
 - Elevate the arm postoperatively for one to three days as necessary
- 4-8 weeks
 - If a greater than 45 degree flexion contracture was present before surgery, use an extension brace at night for 4-8 weeks
- First 3 months
 - The patient must not lift more than one pound (~0.5 kg) during the first three postoperative months
 - Lifting in flexion is not an issue
- After 3 months
 - It is recommended the patient not lift more than 1 kg (2 lbs) on a repetitive basis or more than 2.5 kg (5 lbs) as a single event

Chapter 10: Revision Technique

10.0 Disarticulate the joint

- (Scrub Nurse will prepare the Driver Bit/Removal Driver)
- Remove the two Humeral Screws using the Humeral Screw Removal Driver and Bit (Figure 10.0A)
- Insert the Articulation Extractor between the Humeral Bearing and one of the Ulnar Bearings and gently pry (do not impact) to separate the Articulation (Figure 10.0B)
- Remove the Axle-Pin by using a rongeur to grasp one end of the Axle-Pin and pull it out of the Articulation Assembly with other hand (Figure 10.0C)

Note: An optional Elbow Torque Driver may be used to remove the Humeral Screws





Figure 10.0A



Figure 10.0B



10.1 Ulnar Component Removal

• Insert the Implant Extractor Hook into the eye of the Ulnar Component and remove the ulnar component with the Slide Hammer (Figure 10.1 A&B)







Figure 10.1B

10.2 Humeral Component Removal

4

- Remove the **Humeral Bearing** by inserting the tip of the **Humeral Bearing Extractor** into the hole just below the **Humeral Bearing**
- Gently push (do not impact) the handle down to dislodge the bearing from the base of the **Humeral Component** yoke (Figure 10.2A)




Figure 10.2A

• Attach Humeral Extractor Plate (choose Size 4 or 5/6) by inserting two Humeral Extractor Screws into the holes on the posterior side of the Humeral Component (Figure 10.2B) and lightly tighten with the Small Hex Screwdriver



Note – for the final version of the surgical technique, this image will show the Humeral Extractor Plate being assembled to the Humeral Component, with screws and driver in image.

• Insert the Implant Extractor Hook under the Humeral Stem Extractor Plate and remove with the Slide Hammer (Figure 10.2C)



Figure 10.2C

10.3 Poly Swap

- Place the Humeral Bearing using the Humeral Bearing Placement Tool (Figure 10.3A)
- Fully seat the **Humeral Bearing** by carefully impacting with the **Humeral Bearing Driver** (Figure 10.3B) or alternatively using the **Articulation Inserter** (Figure 10.3C)
 - Confirm no visual gaps between the Humeral Bearing and the Humeral yoke
- Finish procedure using primary technique





Figure 10.3B

This image will show the Articulation Inserter pressing the bearing in place

Figure 10.3C

Chapter 11: Scrub Nurse Instructions



Note: Perform all tasks at the back-table

11.0 Saw Blades

• Attach a .050" thick oscillating saw blade (1.27 mm thick) to power equipment for use during the Humeral Preparation

11.1 Rasp Connections

• Attach the Rasp to the T-Handle (Figure 11.1A)



Figure 11.1A

• Alternatively, if power-rasping is employed, connect one end of the **Reciprocating Saw** Adapter to the drill and the other end to the Single-Use Connector, which then connects to **Rasp** (Figure 11.1D&E)




Figure 11.1D

Attachment to the	
recip saw	

Figure 11.1E

11.2 Humeral Bearing Assembly

- Once Humeral Component size has been determined, place the Humeral Bearing into the Humeral Component using the Humeral Bearing Placement Tool (Figure 11.2A)
 - o The Humeral Bearing does not have a specific anterior or posterior orientation
 - o The Humeral Bearing will not appear to be fully seated at this stage



Note – for the final version of the surgical technique, this image will show the Humeral Bearing being placed into the Humeral Component by the Humeral Placement Tool. Shows handles, no bone in this image (just humeral component).

- Figure 11.2A
- Position the Humeral Bearing Driver against articulation surface of the Humeral Bearing and insert the Humeral Driver Pin (Figure 11.2B)
 - o Keep Driver handle parallel to the Humeral Component rotation axis
 - Allow the loose ring to drop into the pockets of the Humeral yoke to help stabilize the device
- Turn the handle approximately ¹/₄- revolution or about 90-degrees until the **Humeral Bearing** is fully seated
 - o Resistance will be felt, but no audible click will occur
 - **Humeral Bearing** will be fully seated when there are no visual gaps when viewing from the posterior and the anterior sides of the Humeral yoke (Figure 11.2C)
- Remove the Instruments

ŧ



Figure 11.2B



Figure 11.2C

• Alternatively, position the Articulation Inserter against the articulation surface of the Humeral Bearing and insert the peg of the proximal leg of the Articulation Inserter into the hole at the base of the Humeral yoke (Figure 11.2D)




Figure 11.2D


• Squeeze the handles to fully seat the **Humeral Bearing** (Figure 11.2E)

Squeeze handles together

Figure 11.2E

11.3 Ulnar Bearing/Axle-Pin Assembly

- Load an Ulnar Bearing (same box as Humeral Bearing) into one side of the Ulnar Bearing Assembly Tool, and one end of the Axle-Pin into the hole in the opposite jaw of the tool (Figure 11.3)
- Grasp Axle-Pin with your fingers and squeeze the handles together with the other hand until hard resistance is felt




- Load second Ulnar Bearing (DO NOT squeeze the second bearing onto the Axle-Pin)
- Carefully introduce the Ulnar Bearing Assembly Tool, with Ulnar Bearing/Axle and second Ulnar Bearing pre-loaded to the surgeon

11.4 Load Screws

1

- Insert the **Humeral Screws** into the **Humeral Screw Holder** by gently threading them in place using the applied flexible tubing to provide tactility (Figure 11.4)
- Remove and discard tubing from screws





11.5 Prepare the Humeral Screw Installation Driver

- Load the Humeral Screw Installation Driver Bit into the Humeral Screw Installation Driver (Figure 11.5)
- Do not remove the plastic sleeve from the Installation Driver Bits




Figure 11.5

11.6 Prepare the Humeral Screw Removal Driver

8, I

• Load the Humeral Screw Removal Driver Bit into the Humeral Screw Removal Driver (Figure 11.6)

Image shows the connection of the humeral screw removal bit into the humeral screw removal driver

Figure 11.6



Traditional 510(k) Premarket Notification

Sterilization

The Zimmer Nexel Total Elbow is sterilized in the following four groupings: the humeral component, ulnar component, the humeral screws (2 in a package) and the articulation kit (humeral bearing-A, 2 ulnar bearings-B and an axle pin). The humeral component, ulnar component and humeral screws are sterilized using Gamma Irradiation (⁶⁰Co) at a contract sterilizer. The articulation kit is sterilized using Ethylene Oxide Gas processing at a contract sterilizer.

Sterilization Method for Humeral Component, Ulnar Component and Humeral screws

The Zimmer Nexel Total Elbow humeral component, ulnar component and humeral screws will be terminally sterilized by gamma irradiation (⁶⁰Co). The minimum to maximum dose range is 25.0-37.0 kGy. The minimum sterilization dose will be verified using VDmax²⁰ and the gamma irradiation processing and dose mapping were conducted using the following standards:

- ANSI/AAMI/ISO 11137-1:2006, "Sterilization of health care products Radiation Requirements for development, validation and routine control of a sterilization process for medical devices."
- ANSI/AAMI/ISO 11137-2:2006, "Sterilization of health products Radiation Establishing the sterilization dose."
- AAMI TIR 33:2005, "Sterilization of health care products Radiation Substantiation of a selected sterilization dose Method VDmax."

The products will be accepted for release as sterile through the methods stated above to provide a sterility assurance level (SAL) of 10^{-6} or better.

Sterilization Method for the Articulation Kit

The Zimmer Nexel Total Elbow articulation kit, which consists of three polyethylene components (1 humeral bearing-A and 2 ulnar bearings-B) and an axle-pin, will be terminally sterilized by ethylene oxide (EtO) sterilization and sold sterile.

<u>EtO</u> – Sterilization cycles are validated per appropriate requirements of ANSI/AAMI/ISO 11135-1:2007, *Sterilization of health care products – Ethylene Oxide Part 1: Requirements for development, validation, and routine control of a sterilization process for medical devices.* Ethylene Oxide residuals analysis will be conducted to meet the requirements of ANSI/AAMI/ISO 10993-7:2008, *Biological evaluation of medical devices – Part 7: Ethylene oxide sterilization residuals.*

As part of the sterilization validations for EtO sterilization, a residuals analysis dissipation curve test will be conducted on the subject device.



Traditional 510(k) Premarket Notification

The products will be accepted for release as sterile through the methods stated above to provide a sterility assurance level (SAL) of 10^{-6} or better.

Re-sterilization of the polyethylene components (humeral bearing-A and ulnar bearings-B) is not recommended under any circumstance.

The following standards are applicable:

- ANSI/AAMI/ISO 11135:2007, Sterilization of Health care products Ethylene Oxide Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices.
- ANSI/AAMI/ISO 10993-7:2008, Biological Evaluation of Medical Devices Part 7: Ethylene Oxide sterilization residuals
- AAMI TIR 28:2009 Product adoption and process equivalency for Ethylene Oxide sterilization

Shelf Life

The shelf life of the packaging for the Zimmer Nexel Total Elbow humeral component, ulnar component and humeral screws (note: the "humeral screws" are also referred to as the "humeral screw kit") is 10 years. The shelf life of the articulation kit (the humeral bearing-A, humeral bearings-B and the axle pin) is 5 years. This testing was performed according to ASTM F1980 *Standard Guide for Accelerated Aging for Sterile Medical Device Packages* and ANSI/AAMI/ISO 11607:2006 *Packaging for Terminally Sterilized Medical Devices*. The packaging materials demonstrate no degradation over time and the shelf life test data assures sterile packaging integrity.



Biocompatibility

The materials contained in the proposed device along with the standards to which they conform are listed and clinically referenced below:

Material Used	Zimmer Specification (ZES)	Conformance Standard
Zimaloy (Co-Cr-Mo)	2A-74	ASTM F1537
Tivanium (Ti-6Al-4V)	2A-81	ASTM F136
Plasma Spray (Ti-6Al-4V)	2A-51	ASTM F1580
<i>IonGuard</i> ® (nitrogen ion surface enrichment)	2A-65	NA
Vivacit-E TM Vit-E HXPE (Vitamin- E blended highly crosslinked polyethylene)	2B-124	NA

These materials, used in the implants, have been shown to be acceptably compatible with biological tissues through previous clinical usage in orthopedic implant applications:

- Zimaloy is identical to that used in existing, previously cleared components, such as Coonrad-Morrey Total Elbow Ulnar and Humeral Implants manufactured by Zimmer, Inc (K001989, cleared July25,2000)
- *Tivanium* is identical to that used in existing, previously cleared components such as Coonrad-Morrey Total Elbow Ulnar and Humeral Implants manufactured by Zimmer, Inc (K001989, cleared July25,2000)
- Plasma spray is identical to that used in existing, previously cleared components such as Coonrad-Morrey Total Elbow Ulnar Implants manufactured by Zimmer, Inc (K001989, cleared July25,2000)
- *IonGuard*® for Titanium material is identical to that used in existing, previously cleared components such as Natural Knee System manufactured by Zimmer (K853991, cleared Feb21,1986).
- *Vivacit-E*TM Vit-E HXPE is identical to that used in existing, previously cleared components such as *Vivacit E*, Vitamin E Highly Crosslinked Polyethylene Liners manufactured by Zimmer, Inc (K120370, cleared June4, 2012)



Traditional 510(k) Premarket Notification

Software

í

This device is an orthopaedic implant and contains no associated software.




 \mathbf{r}

Traditional 510(k) Premarket Notification

Electromagnetic Compatibility and Electrical Safety

This device does not contain any electronic components; therefore, this section does not apply.



Testing specific to Vivacit- E^{TM} (Vitamin E [α -tocopherol] stabilized, high crosslinked ultra-high molecular weight polyethylene) is contained within the Vitamin E Stabilized HXPE Material Master File authored by Zimmer (FDA Ref: MAF-1868).

Testing specific to plasma spray for the humeral component is contained within the Plasma Sprayed Titanium-6Al-4V Alloy Coating Material Master File (FDA Ref: MAF-1909, Letter of Authorization can be found in Section 21).

Testing specific to plasma spray for the ulnar component is contained within the Coatings Characterization of Orthopedic Titanium 6Al4V Coatings Material Master File (FDA Ref: MAF-652, Letter of Authorization can be found in Section 22).

Consistent with clinically observed failure modes, the Zimmer Nexel Total Elbow test strategy was based on the following keys to clinical success: stem loosening¹, stem fatigue, bearing wear, system durability and modular connection fatigue and resistance to disassociation. Testing methodology was informed by internationally recognized standards, FDA guidance documents, previous testing, biomechanical rationales, peer-reviewed journal articles and predicate device retrieval studies.

Testing was conducted on the proposed device (the Zimmer Nexel Total Elbow) and compared to a legally marketed predicate device (the Zimmer Coonrad/Morrey Total Elbow). The proposed device functioned as intended and met all acceptance criteria, thus demonstrating that the Zimmer Nexel Total Elbow is as safe and effective as the legally marketed Zimmer Coonrad/Morrey Total Elbow.

In addition to conducting the aforementioned device-specific bench testing to substantiate substantial equivalence of the proposed device to the predicate, testing was conducted to evaluate the delamination resistance of Vivacit-E (Vit-E HXPE) as compared to that of conventional GUR 1050 Ultra High Molecular Weight Polyethylene. The conclusions of this testing and the associated proposed marketing claim are presented below, after the device-specific testing summaries.

Testing is summarized below, and the full test reports are included at the end of this section.

¹ Note: Although, it is a recognized mode of failure for total elbow arthroplasty, stem loosening was not the focus of a performance test for the New Zimmer Total Elbow, because the variables (stem geometry, material and surface finish) associated with stem loosening are consistent from the predicate to the proposed device; and, therefore the proposed device (the New Zimmer Total Elbow) is expected to be as safe and effective as the predicate device (the Coonrad/Morrey Total Elbow) in terms of ability to resist stem loosening.



Note: The "Zimmer Nexel Total Elbow" may also be referred to as the "Lunar Elbow" (an internal project name used in lieu of a brand name), "The New Zimmer Total Elbow", "New Zimmer Elbow", "New Elbow", "Zimmer Total Elbow", "Zimmer New Elbow" and the "proposed device". The "Zimmer Coonrad/Morrey Total Elbow" may also be referred to as the "Coonrad/Morrey", the "C/M" and the "predicate device".

Stem Fatigue Testing

Humeral Stem Fatigue Testing Testing Documentation

- Zimmer Research Report ZRR_WA_2554_12, Coonrad/Morrey Elbow Distal Humeral and Proximal Ulnar Implant Stem Fatigue Testing
- Zimmer Research Report ZRR_WA_2589_12, New Zimmer Total Elbow Humeral Implant Distal Stem Fatigue Testing

Cumulative Conclusions of Testing Documentation

- The Zimmer Nexel Total Elbow is as safe and effective as the predicate device (Zimmer Coonrad/Morrey Total Elbow) in terms of humeral stem fatigue strength.
- The fatigue strength of the Zimmer Nexel Total Elbow humeral component stem at (b)(4) load cycles is (b)(4), which is (b)(6 stronger than that of the legally marketed Zimmer Coonrad/Morrey Total Elbow humeral component (b)(4) at a comparable point on the stem.

Background

Due to a lack of fatigue testing standards specific to the humeroulnar joint and because both the proposed device (the Zimmer Nexel Total Elbow) and the predicate device (the Zimmer Coonrad/Morrey Total Elbow) are indicated for use with the same weight-in-hand (5 lbs), the purpose of this set of tests was to ensure that, under worst-case loading conditions, the fatigue strength of the proposed device (Zimmer Nexel Total Elbow) is as safe and effective as that of the predicate device (the Zimmer Coonrad/Morrey Total Elbow).

The humeral component was fatigue tested in the distal region of the stem (the beaded region on the predicate device, and the plasma spray region on the Zimmer Nexel Elbow), because this is consistent with the location of humeral stem fractures observed and reported clinically.

Methods

First, a test was conducted to determine the fatigue strength at (b)(4) load cycles in the beaded region of the predicate device's humeral component.

The Zimmer Nexel Total Elbow humeral component was then fatigue tested in the same manner as the predicate device; and, the fatigue strength of the proposed device was compared to that of the predicate device.



Results

The predicate device humeral component was fatigue tested at various load levels, and the fatigue strength was determined to be 1(b)(4). Example confidence. The full test report can be found in at the end of this section under Zimmer Research Report ZRR_WA_2554_12, Coonrad/Morrey Elbow Distal Humeral and Proximal Ulnar Implant Stem Fatigue Testing.

Five Zimmer Nexel Total Elbow humeral stems achieved (b)(4) load cycles at (b)(4) without fracture. The full test report can be found at the end of this section under Zimmer Research Report ZRR_WA_2589_12, New Zimmer Total Elbow Humeral Implant Distal Stem Fatigue Testing.

Device Update

The following update was made to the humeral component of the Zimmer Nexel Total Elbow device after the conclusion of humeral stem fatigue testing and may be considered theoretically relevant to the outcome of the test, but this update was determined to have no impact on testing results or cumulative conclusions of testing documentation.

Humeral Stem Fatigue Testing, Update 1 (Humeral Component):

The plasma spray pattern has changed slightly, specifically the termination of plasma spray coating has been moved posteriorly to more closely reflect the beaded pattern of the predicate device. See Images 18.1 and 18.2 below.




Image 18.1 Humeral Stem Fatigue Testing, Update 1 (blue line is potting level, gray represents plasma spray region on the **humeral components tested in this study**)




Image 18.2 Humeral Stem Fatgue Testing, Update 2 (blue line is potting level, gray represents plasma spray region on **current Zimmer Nexel Total Elbow** production-intent humeral component)

Rationale for why this update does not change testing results or cumulative conclusions of testing documentation: Because the plasma spray coating reduces the fatigue strength of Tivanium [Zimmer Technical Memo 1222.00, Fatigue Strength of Plasma Sprayed Ti-6AL-4V Alloy Material], moving the plasma spray away from the anterior side (the region that experiences maximum tensile stress) would tend to further reduce the likelihood of a fatigue fracture in the high stress zone; therefore, this update is numerically negligible with respect to the results and conclusions associated with this test.

Internal References

The following research documents, internal to Zimmer, informed this fatigue testing and are provided for reference in their entirety at the end of this section under:

- Zimmer Research Memo ZRM_WA_0292_11, Summary of Clinical Complications of the Coonrad/Morrey Elbow System
- Zimmer Research Report ZRR_WA_2476_11, Biomechanics Rationale for the Loading of Total Elbow Prosthesis During In Vitro Simulations
- Zimmer Technical Memo 1218.00, Coonrad/Morrey Humeral Stem Fatigue
 <u>Testing</u>
- Zimmer Technical Memo 1222.00, Fatigue Strength of Plasma Sprayed Ti-6AL-4V Alloy Material



Ulnar Stem, Mid-Stem Fatigue Testing Testing Documentation

- Zimmer Research Report ZRR_WA_2585_12, Coonrad/Morrey Elbow Ulnar Implant Mid-Stem Fatigue Testing
- Zimmer Research Report ZRR_WA_2615_12, New Zimmer Total Elbow Ulnar Implant Mid-Stem Stem Fatigue Testing

Cumulative Conclusions of Testing Documentation

- The Zimmer Nexel Total Elbow is as safe and effective as the predicate device (Zimmer Coonrad/Morrey Total Elbow) in terms of ulnar stem mid-stem fatigue strength.
- The fatigue strength of the Zimmer Nexel Total Elbow ulnar component stem in the mid-stem region at (b)(4), and load cycles is at least the same as that of the legally marketed Zimmer Coonrad/Morrey Total Elbow ulnar component at a comparable point on the stem (b)(4).

Background

Due to a lack of fatigue testing standards specific to the humeroulnar joint and because both the proposed device (the Zimmer Nexel Total Elbow) and the predicate device (the Zimmer Coonrad/Morrey Total Elbow) are indicated for use with the same weight in hand (5 lbs), the purpose of this set of tests was to ensure that, under worst-case loading conditions, the fatigue strength of the proposed device (Zimmer Nexel Total Elbow) is as safe and effective as that of the predicate device (the Zimmer Coonrad/Morrey Total Elbow).

The ulnar component was subjected to testing along the mid-stem (also referred to as "tapered", "distal", and "non plasma spray") region of the stem to ensure that the Zimmer Nexel Total Elbow ulnar component (which, unlike the predicate device, was subjected to a nitrogen enrichment hardening process) remained as safe and effective as the predicate device in terms of fatigue strength in this mid-stem region.

Methods

First, a test was conducted to determine the fatigue strength in the mid-stem (non plasma spray region) of the predicate device's ulnar component.

The Zimmer Nexel Total Elbow ulnar component was then fatigue tested in the same manner as the predicate device; and, the fatigue strength of the proposed device was compared to that of the predicate device.



186

Results

The predicate device ulnar component was fatigue tested at various load levels, and the fatigue strength in the mid-stem region was determined to be (b)(, as 1 predicate device fractured at (b)(4) achieved (b)(4) load cycles at (b)(without fracture. The full test report can be found at the end of this section under <u>Zimmer</u> <u>Research Report ZRR_WA_2585_12</u>, <u>Coonrad/Morrey Elbow Ulnar Implant Mid-Stem Fatigue Testing</u>.

(b): Zimmer Nexel Total Elbow ulnar stems achieved (b)(4) load cycles in the mid-stem region at a load of (b)(4 without fracture. The full test report can be found at the end of this section under Zimmer Research Report ZRR_WA_2615_12, New Zimmer Total Elbow Ulnar Implant Mid-Stem Stem Fatigue Testing.

Device Updates

The following update was made to the components of the Zimmer Nexel Total Elbow device after the conclusion of ulnar stem mid-stem fatigue testing and may be considered theoretically relevant to the outcome of the test, but these updates were determined to have no impact on testing results or cumulative conclusions of testing documentation.

Ulnar Stem Mid-Stem Fatigue Testing, Update 1 (Ulnar Component):

On the size 4 and size 5 ulnar implants, approximately (b)(4) of material was removed from the anterior side and approximately (b)(4) of material was added to the posterior side of the neck region (the proximal-most aspect of the ulnar stem, closest to the ulnar eye). Similarly, on the size 6 ulnar implant, approximately (b)(4) of material was removed from the anterior side and approximately (b)(4) of material was added to the posterior side of the neck region. This design change was made to increase the overall flexion-extension range of motion to (b)(for all ulnar-humeral implant combinations. See Image 18.3 below.




Image 18.3: Ulnar Stem Mid-Stem Fatigue Testing, Update 1



187

Rationale for why this update does not change testing results or cumulative conclusions of testing documentation: This change reduced the cross-sectional area of the size 4 ulnar component (worst case tested in the ulnar stem mid-stem fatigue test) in the neck region from (b)(4) Because this change does not affect the dimensions at the tested cross-section, it neither alters the moment of inertia nor the moment arm created by the worst-case test orientation; therefore, the stresses at the test location are also not expected to change. Therefore, this update is numerically negligible with respect to the results and conclusions associated with this test.

Internal References

The following research documents, internal to Zimmer, informed this fatigue testing and are provided for reference in their entirety at the end of this section under:

- Zimmer Research Memo ZRM_WA_0292_11, Summary of Clinical Complications of the Coonrad/Morrey Elbow System
- Zimmer Research Report ZRR_WA_2476_11, Biomechanics Rationale for the Loading of Total Elbow Prosthesis During In Vitro Simulations
- Zimmer Technical Memo 1217.00, Coonrad/Morrey Elbow Ulnar Stem Fatigue Testing



Ulnar Stem, Proximal Fatigue Testing Testing Documentation

- Zimmer Research Report ZRR_WA_2554_12, Coonrad/Morrey Elbow Distal Humeral and Proximal Ulnar Implant Stem Fatigue Testing
- Zimmer Research Report ZRR_WA_2626_12, New Zimmer Total Elbow Ulnar Implant Proximal Stem Fatigue Testing

Cumulative Conclusions of Testing Documentation

- The Zimmer Nexel Total Elbow is as safe and effective as the predicate device (Zimmer Coonrad/Morrey Total Elbow) in terms of ulnar component stem fatigue strength in the plasma spray (proximal stem) region.
- The fatigue strength of the Zimmer Nexel Total Elbow ulnar component stem in the plasma spray (proximal stem) region at (b)(4) load cycles is at least the same as that of the legally marketed Zimmer Coonrad/Morrey Total Elbow ulnar component at a comparable point on the stem (b)(4), and thus meets the performance requirement for substantial equivalence for the ulnar component stem in the plasma spray region.
- The fatigue strength of the Zimmer Nexel Total Elbow ulnar component stem in the plasma spray (proximal stem) region at (b)(4) load cycles is (b)(4 (a load equivalent to (b)(4) weight-in-hand), which exceeds the performance requirement for substantial equivalence for the ulnar component stem in the plasma spray region.

Background

Due to a lack of fatigue testing standards specific to the humeroulnar joint and because both the proposed device (the Zimmer Nexel Total Elbow) and the predicate device (the Zimmer Coonrad/Morrey Total Elbow) are indicated for use with the same weight in hand (5 lbs), the purpose of this set of tests was to ensure that, under worst-case loading conditions, the fatigue strength of the proposed device (Zimmer Nexel Total Elbow) is as safe and effective as that of the predicate device (the Zimmer Coonrad/Morrey Total Elbow).

The ulnar component was fatigue tested in the plasma spray region of the stem, because this is consistent with the location of ulnar stem fractures observed and reported clinically.

Methods

First, a test was conducted to determine the fatigue strength of the predicate device ulnar component in the plasma spray (proximal stem) region.

The Zimmer Nexel Total Elbow ulnar component was then fatigue tested in the same manner as the predicate device; and, the fatigue strength of the proposed device was compared to that of the predicate device.



Results


The predicate device ulnar component was fatigue tested at various load levels, and the fatigue strength in the plasma spray (proximal stem) region was determined to be b) as 1 predicate device fractured at (b)(1 fractured at (b)(1 fractured at and (b)(4) load cycles at (b)(4 without fracture. The full test report can be found at the end of this section under <u>Zimmer Research Report</u> <u>ZRR_WA_2554_12, Coonrad/Morrey Elbow Distal Humeral and Proximal Ulnar</u> Implant Stem Fatigue Testing.

(b) (Zimmer Nexel Total Elbow ulnar stems achieved (b)(4) load cycles in the plasma spray (proximal stem) region at a load of (b)(4) without fracture, meeting the performance requirement. In addition, five Zimmer Nexel Total Elbow ulnar stems achieved (b)(4) load cycles in the plasma spray (proximal stem) region at a load of (b)(4) (a load equivalent to (b)(4) los weight-in-hand) without fracture. The full test report can be found at the end of this section under Zimmer Research Report ZRR WA 2626 12. New Zimmer Total Elbow Ulnar Implant Proximal Stem Fatigue Testing.

Device Updates

The following update was made to the components of the Zimmer Nexel Total Elbow device after the conclusion of ulnar stem proximal fatigue testing and may be considered theoretically relevant to the outcome of the test, but these updates were determined to have no impact on testing results or cumulative conclusions of testing documentation.

Ulnar Stem Proximal Fatigue Testing, Update 1 (Ulnar Component): On the size 4 and size 5 ulnar implants, approximately (b)(4), of material was removed from the anterior side and approximately (b)(4), of material was added to the posterior side of the neck region (the proximal-most aspect of the ulnar stem, closest to the ulnar eye). Similarly, on the size 6 ulnar implant, approximately (b)(4) of material was removed from the anterior side and approximately (b)(4) of material was added to the posterior side of the neck region. This design change was made to increase the overall flexion-extension range of motion to (b)(for all ulnar-humeral implant combinations. See Image 18.4 below.



Rationale for why this update does not change testing results or cumulative conclusions of testing documentation: This change reduced the cross-sectional area of the size 4 ulnar component (worst case tested in the ulnar stem proximal fatigue test) in the neck region from (b)(4). Because this change does not affect the dimensions at the tested cross-section (in the plasma spray region), it neither alters the moment of inertia nor the moment arm created by the worst-case test orientation; therefore, the stresses at the test location (in the plasma spray region) is also not expected to change.

This rationale for why this device update does not create a new worst case in the ulnar stem proximal fatigue testing is further substantiated by an additional finite element analysis performed (ZRR_WA_2648_12), which concluded that the peak tensile stresses (during similar clinically-relevant worst-case loading conditions) in the neck region of the production-intent Nexel ulnar component remain below the experimentally demonstrated material fatigue strength of wrought Tivanium. Therefore, this update is numerically negligible with respect to the results and conclusions associated with this test.

Internal References

The following research documents, internal to Zimmer, informed this fatigue testing and are provided for reference in their entirety at the end of this section under:

- Zimmer Research Memo ZRM_WA_0292_11, Summary of Clinical Complications of the Coonrad/Morrey Elbow System
- Zimmer Research Report ZRR_WA_2476_11, Biomechanics Rationale for the Loading of Total Elbow Prosthesis During In Vitro Simulations
- Zimmer Technical Memo 1217.00, Coonrad/Morrey Elbow Ulnar Stem Fatigue Testing
- Zimmer Research Report ZRR_WA_2648_12, Finite Element Analysis of Nexel Elbow Ulnar Implant Neck Region



Wear Testing

Testing Documentation

- Zimmer Research Report ZRR_WA_2407_11, Coonrad/Morrey Wear Test Method
- Zimmer Research Report ZRR_WA_2552_12, Wear Testing of the New Zimmer Total Elbow Prosthesis

Cumulative Conclusions of Testing Documentation

- The Zimmer Nexel Total Elbow met or exceeded all acceptance criteria and is at least as safe and effective as the predicate device in terms of wear performance.
- The Zimmer Nexel Total Elbow mean gravimetric wear rate (b)(4) at (b)(4) load cycles) was (b) less than that of the predicate.

Background

Due to a lack of wear testing standards specific to the humeroulnar joint, a wear test was developed by Zimmer that produces results on the Zimmer Coonrad/Morrey Total Elbow polyethylene bearings similar to those observed on clinically retrieved Zimmer Coonrad/Morrey Total Elbow polyethylene bearings (the full test report can be found at the end of this section under <u>Zimmer Research Report</u> <u>ZRR WA 2407 11, Coonrad/Morrey Wear Test Method</u>).

The same wear test was then conducted using the proposed (Zimmer Nexel Total Elbow) device in order to ensure that the proposed device is equal to or better than the predicate device (the Zimmer Coonrad/Morrey Total Elbow) in terms of gravimetric wear results. The full test report can be found at the end of this section under <u>Zimmer</u> <u>Research Report ZRR_WA_2552_12</u>, <u>Wear Testing of the New Zimmer Total Elbow</u> <u>Prosthesis</u>.

Methods

(b)(4) proposed device (Zimmer Nexel Total Elbow) articular assemblies were tested at (b)(4) with a loading profile shown to produce clinically relevant results for the predicate device (the Zimmer Coonrad/Morrey Total Elbow). This testing on the proposed device was conducted in four intervals and the gravimetric wear rate determined at (b)(4) cycles. In addition, the test lubricant was salvaged from (b)(4) randomly selected stations at the (b)(4) cycle milestones and the polyethylene debris isolated and compared to the debris obtained from the same test conducted on the predicate device.

The Zimmer Nexel Total Elbow was required to complete this test regime without component fracture, device disassembly, or evidence of metal-on-metal articulation. Additionally, the Zimmer Nexel Total Elbow was required to complete this test with a



mean gravimetric wear rate equal to or less than that of the predicate device (b)((b)(4)). Finally, all humeral screws were required to have measureable torque upon removal at the conclusion of testing.

Results

Al (b)(4) Zimmer Nexel Total Elbow test specimens completed the testing protocol without incidences of component fracture, disassembly or evidence of metal-on-metal articulation and with isolated poly debris similar in shape (aspect ratio, AR) and size (equivalent circle diameter, ECD) as compared to the poly debris isolated from the predicate device wear test. The Zimmer Nexel Total Elbow mean gravimetric wear rate (b)(4) load cycles) was(b)(less than that of the predicate device (b)(4) load cycles). All humeral screws had measureable torque upon removal at the conclusion of testing.

Device Updates

The following updates were made to the components of the Zimmer Nexel Total Elbow device after the conclusion of wear testing and may be considered theoretically relevant to the outcome of the test, but these updates were determined to have no impact on testing results or cumulative conclusions of testing documentation.

Wear Testing, Update 1 (Ulnar Component):

The width of the flat surface on the inner diameter of the ulnar eye increased slightly from (b)(4) (causing an overall inner diameter increase of the ulnar eye by (b)(4) in order to improve load distribution and ensure the allowable minimum pistoning (proximodistal and anteroposterior translational) motion of the proposed device is no less than that of the predicate device. See Image 18.5 below.




Image 18.5: Wear Testing, Update 1

Rationale for why this update does not change testing results or cumulative conclusions of testing documentation: The surfaces affected by this update did not articulate against the ulnar bearings-B components during this wear test as a result of



the fixed varus/valgus position of the ulnar component with respect to the poly and humeral components dictated by the clinically-relevant worst-case test method; therefore, this update is numerically negligible with respect to the results and conclusions associated with this test.

Wear Testing, Update 2 (Ulnar Component):

The mediolateral radial dimension of the ulnar eye was increased by (b)(4) (causing a decrease of the anteroposterior outer diameter of the ulnar eye by (b)(4) to improve load distribution. See Image 18.6 below.





Image 18.6: Wear Testing, Update 2

Rationale for why this update does not change testing results or cumulative conclusions of testing documentation: This update is numerically negligible with respect to the results and conclusions associated with this test.

Wear Testing, Update 3 (Ulnar Bearings-B):

The compressible semi-circle design feature was originally only on one of the ulnar bearings (ulnar bearing-B). This compressible semi-circle was then added to the second ulnar bearing (originally called ulnar bearing-C) in order to make the two ulnar bearings identical and avoid the potential use error of choosing the wrong of two nearly identical components. See Image 18.7 below.





Rationale for why this update does not change testing results or cumulative conclusions of testing documentation: This update only very slightly improved the overall press-fit between apposing ulnar bearings-B and between ulnar bearings-B and the humeral yoke, making this update numerically negligible with respect to the results and conclusions associated with this test.

Wear Testing, Update 4 (Ulnar Bearings-B):

The distance between the major external face and the minor internal face was decreased by (b)(4) to ensure non-interference between the broad faces of apposing ulnar bearings-B. See Image 18.8 below.




Image 18.8: Wear Testing, Update 4



Rationale for why this update does not change testing results or cumulative conclusions of testing documentation: This update does not affect any articulating surface, and the press-fit of tested assemblies is representative of the current production-intent device.

Wear Testing, Update 5 (Ulnar Bearings-B):

The overall diametric size of the central (shoulder) articulation region of ulnar bearings-B was reduced by (b)(4) n order to ensure the allowable minimum pistoning (proximodistal and anteroposterior translational) motion of the proposed device is no less than that of the predicate device. See Image 18.9 below.




Image 18.9: Wear Testing, Update 5



Rationale for why this update does not change testing results or cumulative conclusions of testing documentation: This slight reduction in the nominal cross-sectional size of the central region of the ulnar bearings-B is numerically negligible with respect to the results and conclusions associated with this test.


Wear Testing, Update 6 (Ulnar Component):

The nitrogen-enrichment process used to increase the hardness of the Zimmer Nexel Total Elbow ulnar component articulating surfaces (the ulnar eye) was changed from Ti-Nidium® (furnace nitriding) to IonGuard® (nitrogen ion implantation).

Rationale for why this update does not change testing results or cumulative conclusions of testing documentation: This update is a process change resulting in the same surface hardness and surface roughness specifications, and is therefore numerically negligible with respect to the results and conclusions associated with this test. This rationale is further substantiated by the pin-on-flat wear test conducted (ZRR_WA_2646_12), which concluded that there is no statistical difference between the average weight loss of a VitE HXPE polymer pin sample coupled with a Ti-Nidium® treated disc and a VitE HXPE polymer pin sample coupled with an IonGuard® treated disc.

Wear Testing, Update 7 (Ulnar Component):

On the size 4 and size 5 ulnar implants, approximately (b)(4) of material was removed from the anterior side and approximately (b)(4) of material was added to the posterior side of the neck region (the proximal-most aspect of the ulnar stem, closest to the ulnar eye). Similarly, on the size 6 ulnar implant, approximately (b)(4)of material was removed from the anterior side and approximately (b)(4) of material was added to the posterior side of the neck region. This design change was made to increase the overall flexion-extension range of motion to (b)(-) or all ulnar-humeral implant combinations. See Image 18.10 below.






Image 18.10: Wear Testing, Update 7



Internal References

The following research documents, internal to Zimmer, informed this wear testing and are provided for reference in their entirety at the end of this section under:

- Zimmer Research Memo ZRM_WA_0292_11, Summary of Clinical Complications of the Coonrad/Morrey Elbow System
- Zimmer Research Report ZRR_WA_2476_11, Biomechanics Rationale for the Loading of Total Elbow Prosthesis During In Vitro Simulations
- <u>Zimmer Research Memo ZRM_WA_0193_09</u>, Development of the <u>Coonrad/Morrey Elbow Durability Test</u>
- Zimmer Research Memo ZRM_WA_0179_09, Analysis of Retrieved Coonrad/Morrey Elbow Total Arthroplasty Components
- Zimmer Research Report ZRR_WA_1988_09, Literature Review Summary of Elbow Joint Mechanics and Total Elbow Arthroplasty Outcomes
- Zimmer Research Report ZRR_WA_2338_11, Vitamin E Highly Crosslinked Polyethylene Shelf Life Study
- Zimmer Research Report ZRR_WA_2646_12, Pin-on-Flat Wear Comparison
 of the New Zimmer Total Elbow and the Coonrad-Morrey Articulation
 <u>Couples</u>



Durability Testing

Testing Documentation

- Zimmer Research Memo ZRM_WA_0300_11, Total Elbow Durability Test
 Methodology Development
- Zimmer Research Report ZRR_WA_2542_12, Durability Testing of <u>Coonrad/Morrey Total Elbow and New Zimmer Total Elbow</u>

Cumulative Conclusions of Testing Documentation

- The Zimmer Nexel Total Elbow device achieved (b)(4)
 run-outs at the equivalent of (b)(4)
 weight-in-hand and (b)(4)
 varus-valgus moment (b)(4)
- The Zimmer Nexel Total Elbow is at least as safe and effective as the predicate device (Zimmer Coonrad/Morrey Total Elbow) in terms of articular construct durability.
- The Zimmer Nexel Total Elbow is more durable than the predicate device (Zimmer Coonrad/Morrey Total Elbow).

Background

In addition to wear testing the proposed device (the Zimmer Nexel Total Elbow) with a relatively low weight-in-hand and relatively high load frequency (replicating normal activities of daily living), Zimmer conducted an additional (durability) elbow prosthesis performance test to assess resistance to gross mechanical failure under strenuous use (relatively high weights-in-hand at a relatively low frequency of occurrence) conditions.

Due to a lack of durability testing standards specific to the humeroulnar joint, a test was developed by Zimmer, using the predicate device (the Zimmer Coonrad/Morrey Total Elbow), that reproduces clinical complications (such as prosthesis disassembly, severe ulnar and humeral bushing/bearing wear resulting in metal-on-metal articulation and burnishing patterns on the snap/axle-pin components) associated with durability-type failure modes observed on clinically failed and retrieved predicate device assemblies. This test method development can be found in its entirety at the end of this section under Zimmer Research Memo ZRM_WA_0300_11, Total Elbow Durability Test Methodology Development.

This durability test was conducted using both the predicate device and the Zimmer Nexel Total Elbow in order to ensure that the proposed device is equal to or better than the predicate device in terms of durability performance. The full test report can be found in at the end of this section under <u>Zimmer Research Report</u> <u>ZRR_WA_2542_12</u>, <u>Durability Testing of Coonrad/Morrey Total Elbow and New</u> <u>Zimmer Total Elbow</u>.



Methods

Worst-case size predicate device test specimens were tested under physiologically relevant loading conditions at various load levels at (b)(4) de-ionized water in order to obtain (b)(4) run-outs to (b)(4) cycles at a given load level. The run out load for the predicate device was found to be equivalent to (b)(4) weight-in-hand (b)(4) peak joint reaction force [JRF] and (b)(4) varus-valgus moment).

(b)(4) worst-case size Zimmer Nexel Total Elbow test specimens were then tested under the same conditions at the equivalent of (b)(4) weight-in-hand (b)(4) peak JRF and (b)(4) varus-valgus moment), a load level (b)(higher than the load level at which the predicate device ran out. Additionally, all humeral screws were required to have measureable torque upon removal at the conclusion of testing.

Results

The Zimmer Nexel Total Elbow Humeral Screws did not fracture or exhibit stripping during insertion to the specified minimum and maximum torque values (b)(4) run-outs were obtained at an equivalent o (b)(4) weight-in-hand (a load level one (b)(5) step higher than the load level at which the C/M Elbow ran out), and there was no evidence of the following durability-type failures prior to or after reaching (b)(4) cycles: metal-on-metal articulation, prosthesis assembly linkage dissociation, fracture of any component, nor protrusion of the Humeral screws (screw back-out) from their holes in the Humeral Test Specimens. Additionally, after completion of (b)(4) cycles, the test specimens were able to be completely disassembled and the Humeral Screws had measurable torque during removal.

Device Updates

The following updates were made to the components of the Zimmer Nexel Total Elbow after the conclusion of durability testing and may be considered theoretically relevant to the outcome of the test, but these updates were determined to have no impact on testing results or cumulative conclusions of testing documentation.

Durability Testing, Update 1 (Ulnar Component):

The width of the flat surface on the inner diameter of the ulnar eye increased slightly from (b)(4) (causing an overall inner diameter increase of the ulnar eye by (b)(4) to improve load distribution and ensure the allowable minimum pistoning (proximodistal and anteroposterior translational) motion of the proposed device is no less than that of the predicate device. See Image 18.11 below.

Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017




Image 18.11: Durability Testing, Update 1

Rationale for why this update does not change testing results or cumulative conclusions of testing documentation: The slight improvement in load distribution associated with this update is numerically negligible with respect to the results and conclusions associated with this test.

Durability Testing, Update 2 (Ulnar Component):

The mediolateral radial dimension of the ulnar eye was increased by (b) " (causing a decrease of the anteroposterior outer diameter of the ulnar eye by (4) ") to improve load distribution. See Image 18.12 below




Image 18.12: Durability Testing, Update 2

Rationale for why this update does not change testing results or cumulative conclusions of testing documentation: The slight improvement in load distribution associated with this update is numerically negligible with respect to the results and conclusions associated with this test.

Traditional 510(k) Premarket Notification



Durability Testing, Update 3 (Ulnar Bearings-B):

The overall diametric size of the central articulation region of ulnar bearings-B was reduced by (b)(4), in order to ensure the allowable minimum pistoning (proximodistal and anteroposterior translational) motion of the proposed device is no less than that of the predicate device. See Image 18.13 below.




Image 18.13: Durability Testing, Update 3

Rationale for why this update does not change testing results or cumulative conclusions of testing documentation: This slight reduction in the nominal cross-sectional size of the central region of the ulnar bearings-B is numerically negligible with respect to the results and conclusions associated with this test.



Durability Testing, Update 4 (Ulnar Bearings-B):

Slightly modified compound radial geometry on the outer surface (shoulder) of the ulnar bearings-B to ensure proper seating of the bearings within the humeral yoke in all possible tolerance stack-up conditions. See Image 18.14 below.




Image 18.14: Durability Testing, Update 4

Rationale for why this update does not change testing results or cumulative conclusions of testing documentation: This slight geometrical change is numerically negligible with respect to the results and conclusions associated with this test.

Durability Testing, Update 5 (Ulnar Component):

The nitrogen-enrichment process used to increase the hardness of the Zimmer Nexel Total Elbow ulnar component articulating surfaces (the ulnar eye) was changed from Ti-Nidium® (furnace nitriding) to IonGuard® (nitrogen ion implantation).

Rationale for why this update does not change testing results or cumulative conclusions of testing documentation: This update is a process change resulting in the same surface hardness and surface roughness specifications, and is therefore numerically negligible with respect to the results and conclusions associated with this test.



Durability Testing, Update 6 (Ulnar Component)

On the size 4 and size 5 ulnar implants, approximately (b)(4) of material was removed from the anterior side and approximatel (b)(4) of material was added to the posterior side of the neck region (the proximal-most aspect of the ulnar stem, closest to the ulnar eye). Similarly, on the size 6 ulnar implant, approximately (b)(4) of material was removed from the anterior side and approximately (b)(4) of material was added to the posterior side of the neck region. This design change was made to increase the overall flexion-extension range of motion to (b)(for all ulnar-humeral implant combinations. See Image 18.15 below.




Image 18.15: Durability Testing, Update 6

Rationale for why this update does not change testing results or cumulative conclusions of testing documentation: This change reduced the lateral articulation contact area from (b)(4) to (b)(4) and reduced the medial articulation contact area from (b)(4) to (b)(4) and reduced the medial articulation contact area from (b)(4) to (b)(4) these two areas articulate on the major internal faces of the ulnar bearings (bearings B and C, as defined in Wear Testing Update 3, above) during Varus/Valgus movements of the ulnar relative to the humerus. Given the fact that this change in the outermost contact surfaces of the ulnar articulation is slight, compounded with the fact that the central and medial/lateral contoured surfaces of the ulnar eye carry and distribute the vast majority of the loading onto the apposing articulating surfaces of the ulnar bearings B (as opposed to the outermost contact surfaces affected by this change), this update is expected to be numerically negligible with respect to the results and conclusions associated with this test.



Internal References

The following research documents, internal to Zimmer, informed this durability testing and are provided for reference in their entirety at the end of this section under:

- <u>Zimmer Research Memo ZRM_WA_0292_11, Summary of Clinical</u> <u>Complications of the Coonrad/Morrey Elbow System</u>
- Zimmer Research Report ZRR_WA_2476_11, Biomechanics Rationale for the Loading of Total Elbow Prosthesis During In Vitro Simulations
- Zimmer Research Memo ZRM_WA_0179_09, Analysis of Retrieved <u>Coonrad/Morrey Elbow Total Arthroplasty Components</u>
- Zimmer Research Report ZRR_WA_1988_09, Literature Review Summary of Elbow Joint Mechanics and Total Elbow Arthroplasty Outcomes
- Zimmer Research Report ZRR_WI_2441_11, Finite Element Analysis of the Coonrad-Morrey Total Elbow System
- Zimmer Research Report ZRR_WA_2338_11, Vitamin E Highly Crosslinked Polyethylene Shelf Life Study



Modular Connection Fatigue Testing

Testing Documentation

Zimmer Research Report ZRR_WA_2598_12, New Zimmer Total Elbow Prosthesis Modular Connection Fatigue Testing

Cumulative Conclusions of Testing Documentation

- The Zimmer Nexel Total Elbow is safe and effective in terms of nonarticulating, mechanically locked, modular connection fatigue.
- The median fatigue strength of the non-articulating, mechanically locked, modular Zimmer Nexel Total Elbow implant components was no less than
 (b)(4).

Background

Modular junctions have an inherent risk of adverse events that have the potential to compromise the system fatigue strength resulting from loosening, micromotion, fretting and galvanic corrosion due to contact of dissimilar materials.

Because the predicate device (the Zimmer Coonrad/Morrey Total Elbow) does not have non-articulating, mechanically-locked, modular implant components comparable to those of the proposed device (the Zimmer Nexel Total Elbow), side-by-side testing comparing the predicate device to the proposed device was irrelevant. Instead, testing recommended in the Draft FDA Guidance Document dated 1 May 1995 and entitled *Guidance Document for Testing Non-Articulating, "Mechanically Locked", Modular Implant Components* was conducted to provide reasonable assurance of safety and effectiveness of the Zimmer Nexel Total Elbow in terms of modular connection fatigue.

The loading level for this test was derived from the median proximal (plasma spray region) <u>fatigue strength</u> of the predicate device ulnar stem.

Methode (b)(4)

(b)(4) worst-case size Zimmer Nexel Total Elbow non-articulating, mechanically locked, modular connections were subjected to cyclic loading under worst-case physiologic conditions for (b)(4) Testing was conducted at a frequency of (b)(4) a Ringer's solution maintained at (b)(4) and a pH (b)(4). In addition to reporting mean system mass loss of the modular connection, occurrence (or non-occurrence) of the following three key failure modes commonly associated with modular connections were recorded: component fracture, component disassociation and lack of measurable screw torque (i.e. screw back-out).

Results

The median fatigue strength of the non-articulating, mechanically locked, modular Zimmer Nexel Total Elbow implant components was no less than **45**/40. No components fractured or disassociated prior to the end of the test; and, all humeral



screws had measureable torque upon removal at the conclusion of testing. The mean system mass loss of the loaded and soak-control Zimmer Nexel Total Elbow specimens was (b)(4) get the provided and soak respectively.

Internal References

The following research documents, internal to Zimmer, informed this modular connection fatigue testing and are provided for reference in their entirety at the end of this section under:

- Zimmer Research Memo ZRM_WA_0292_11, Summary of Clinical Complications of the Coonrad/Morrey Elbow System
- Zimmer Research Report ZRR_WA_2476_11, Biomechanics Rationale for the Loading of Total Elbow Prosthesis During In Vitro Simulations
- Zimmer Research Report ZRR_WA_1988_09, Literature Review Summary of Elbow Joint Mechanics and Total Elbow Arthroplasty Outcomes
- Zimmer Research Report ZRR_WA_2481_11, Determination of the Functional Torque Range of the New Zimmer Total Elbow Humeral Screw and Hex Driver
- <u>Memo to DHF File 207-014, Rationale for Determination of Torque</u> <u>Specification and Torque Range Requirement for the Mechanical Fastening</u> <u>Used in Zimmer New Total Elbow Prosthesis</u>
- Zimmer Research Report ZRR_WA_2252_10, Evaluation of the AS/BF Taper Adaptor Component Using Electrochemical Test Methods




Bench Testing in support of Vivacit-E Claim

Proposed Claim

Vivacit-E Vitamin E Highly Crosslinked Polyethylene exhibits improved delamination resistance compared to conventional polyethylene under aged conditions (per ASTM F2003). The results of *in vitro* delamination tests have not been shown to correlate with clinical delamination mechanisms.

Evidence to Support Claim



Records processed under FOIA request #2016-4653; Released by CDRH on 07/16/2017



The above claim is supported by the following references:

[1]. <u>ZRR_WA_2580_12</u>, Evaluation of the Delamination Resistance of *Vivacit-E*TM Ultra High Molecular Weight Polyethylene.

[2]. Goldberg SH, Urban RM, Jacobs JJ, et al. Modes of Wear After Semiconstrained Total Elbow Arthroplasty. J Bone Joint Surg Am. 2008; 90:609-619.

[3]. Carol J. Bell, Peter S. Walker, Melanie R. Abeysundera et al. Effect of oxidation on delamination of ultrahigh-molecular-weight polyethylene tibial components. J Arthroplasty. 1998; 13, Issue 3: 280-290

[4]. <u>ZRR_WA_2409_11</u>, Extended Aging of Vitamin E Highly Crosslinked Polyethylene (Master File volume 2, page 394).

[5].<u>ZRR_WA_2399_11_</u>REV1, Effect of Prolonged Aging on the Wear Performance of Vitamin E IT Liner. (Testing from predicate *Vivacit-E* Vitamin E Highly Crosslinked Polyethylene Liners 510(k) submission K120370.)

[6]. <u>TM1267</u>, Delamination Sub-Surface Fatigue Wear Test of Electron-Beam Crosslinked Polyethylene (Testing from *Prolong* Delamination Claim Special 510(k) submission K013991).

[7].D'Lima DD, Steklov N, Fregly B, et al. In Vivo Contact Stresses During Activities of Daily Living After Knee Arthroplasty. J Orthop Res 26:1549-1555 (2008).

[8]. ZRR_WI_1222_12, Finite Element Analyses of the new Zimmer Total Elbow.

Page 1 of 15



Zimmer Research Report

