



**510(K) SUMMARY**

**MAR 24 2013**

***FiberLase Endure CO2 Fiber***  
**510(k) Number K 130164**

**Applicant's Name:** Lumenis Ltd.

13 Hayetzira Street  
Yokneam Industrial Park  
Yokneam 20692 Israel  
Tel. (972)4-959-9000  
Fax: (972)4-959-9050

**Contact Person:** Yoram Levy, Qsite

31 Haavoda St.  
Binyamina, Israel 30500  
Tel (972)4-638-8837  
Fax (972)4-638-0510  
Yoram@qsitemed.com

**Trade Name:** *FiberLase Endure CO2 Fiber*

**Device Type:** *CO2 Laser fiber*

**Preparation Date:** January 20, 2013

**Classification:** **Regulatory Name:** Laser surgical instrument for use in general and plastic surgery and in dermatology

**Product Code:** GEX

**Regulation No:** 21 CFR 878.4810

**Class:** II

**Classification Panel:** General & Plastic Surgery

**Device Description:**



The *FiberLase Endure CO2 Fiber* is a hollow, semi rigid, light-conducting delivery fiber, 2 meter length, designed to transmit laser energy from the CO2 laser system to the treatment site. The fiber also transmits a low power red diode or helium neon laser aiming beam to assist in targeting the tissue to be treated.

**Intended Use Statement:**

The *FiberLase Endure CO2 Fiber* is intended for use in surgical procedures requiring ablation, vaporization, excision, incision and coagulation of soft tissue. The *FiberLase Endure CO2 Fiber* is indicated for use in open surgical procedures such as ENT surgery and laparoscopy and endoscopic procedures.

**Predicate Devices:** Substantial equivalence to the following predicate device is claimed:

Device Name	510k No	Date of Clearance
FiberLase CO2 Laser WaveGuide	K100384	April 12 , 2010

**Performance Standards**

*FiberLase Endure CO2 Fiber* was tested and complies with the following standards:

- ANSI/AAMI/ISO 11135-1:2007 Sterilization of health care products — Ethylene oxide
- AAMI TIR30:2001 Product adoption and process equivalency for ethylene oxide sterilization
- ISO 14971-1:2007 Risk management for medical devices
- IEC 60601-2-22 ed3.0:2007 – Medical Electrical Equipment - Part 2-22: Particular Requirements for Basic Safety and Essential Performance of Surgical, Cosmetic, Therapeutic and Diagnostic Laser Equipment.
- IEC 60825-1:2007 – Safety of Laser Products - Part 1: Equipment Classification, Requirements and User's Guide
- ANSI/AAMI/ISO 17665-1:2006 Sterilization of health care products — Moist Heat.



A detailed description follows in **Section 14**.

### **Performance Testing**

Performance testing demonstrated that the *FiberLase Endure CO2 Fiber* is as safe and effective as the cleared predicate device.

### **Comparison with the Predicate Devices**

The *FiberLase Endure CO2 Fiber* is a modification to its predicate device, the FDA-cleared Lumenis Fiberlase CO2 laser WaveGuide (K100384).

The intended use of the *FiberLase Endure CO2 Fiber* is identical to the intended use of its predicate.

Both the *FiberLase Endure CO2 Fiber* and the Lumenis Fiberlase CO2 laser WaveGuide systems are fibers that transmit laser energy from the laser system to the treatment site. Both devices are comprised of a laser connector and a 2 meter long delivery fiber.

The structures, the materials and the dimensions of the *FiberLase Endure CO2 Fiber* are identical to the cleared Lumenis Fiberlase CO2 laser WaveGuide fiber.

The minor difference between the *FiberLase Endure CO2 Fiber* System and its predicate device is an addition of an Autoclave sterilization method for reprocessing the fiber instead of single use ETO. The minor differences do not raise any new questions of safety or efficacy. Moreover, performance testing demonstrated that the *FiberLase Endure CO2 Fiber* is as safe and effective as the predicate device. Thus, the *FiberLase Endure CO2 Fiber* is substantially equivalent to Lumenis FiberLase CO2 Laser WaveGuide (K100384).



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

Lumenis, Limited  
% Qsite  
Mr. Yoram Levy  
General Manager  
31 Haavoda Street  
Binyamina, Israel 30500

March 24, 2012

Re: K130164

Trade/Device Name: FiberLase Endure CO2 Fiber  
Regulation Number: 21 CFR 878.4810  
Regulation Name: Laser surgical instrument for use in general and plastic surgery and in dermatology  
Regulatory Class: Class II  
Product Code: GEX  
Dated: March 06, 2013  
Received: March 11, 2013

Dear Mr. Levy:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

# LUMENIS

## INDICATIONS FOR USE STATEMENT

510(k) Number (if known): K130164

Device Name: *FiberLase Endure C02 Fiber*

Indications for Use: The *FiberLase Endure C02 Fiber* is intended for use in surgical procedures requiring ablation, vaporization, excision, incision and coagulation of soft tissue. *FiberLase Endure C02 Fiber* is indicated for use in open surgical procedures such as ENT surgery and laparoscopy and endoscopic procedures.


Prescription Use:  (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use:  (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE -CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Neil R Ogden   
2013.03.22 16:12:33-04'00'

(Division Sign-Off) for MXM  
Division of Surgical Devices  
510(k) Number K130164

1 of 1

*FiberLase ENDURE C02 Fiber-5/0k Notification*

K130164

v.1



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

Lumenis, Limited  
% Qsite  
Mr. Yoram Levy  
General Manager  
31 Haavoda Street  
Binyamina, Israel 30500

March 24, 2012

Re: K130164

Trade/Device Name: FiberLase Endure CO2 Fiber  
Regulation Number: 21 CFR 878.4810  
Regulation Name: Laser surgical instrument for use in general and plastic surgery and in dermatology  
Regulatory Class: Class II  
Product Code: GEX  
Dated: March 06, 2013  
Received: March 11, 2013

Dear Mr. Levy:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

1

Page 2 – Mr. Yoram Levy

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/cdrh/mdr/> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,  
FOR

**Peter D. Rumm -S**

Mark N. Melkerson  
Acting Director  
Division of Surgical Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

Page 3 – Mr. Yoram Levy

Concurrence & Template History Page  
 [THIS PAGE IS INCLUDED IN IMAGE COPY ONLY]


**Full Submission Number:**

For Office of Compliance Contact Information:

[http://insideportlets.fda.gov:9010/portal/page?\\_pageid=197.415881&\\_dad=portal&\\_schema=PORTAL&org=318](http://insideportlets.fda.gov:9010/portal/page?_pageid=197.415881&_dad=portal&_schema=PORTAL&org=318)

For Office of Surveillance and Biometrics Contact Information:

[http://insideportlets.fda.gov:9010/portal/page?\\_pageid=197.415881&\\_dad=portal&\\_schema=PORTAL&org=423](http://insideportlets.fda.gov:9010/portal/page?_pageid=197.415881&_dad=portal&_schema=PORTAL&org=423)

<b>Digital Signature Concurrence Table</b>	
Reviewer Sign-Off	Xin (Sofia) Tan, Ph.D.
Branch Chief Sign-Off	Neil R.P. Ogden, MS
Division Sign-Off	Peter D. Rumm -S  2013.03.25 22:15:12-04:00

f/t:XNT:kdm:3/25/13

**Template Name:** K1(A) – SE after 1996

**Template History:**

Date of Update	By	Description of Update
7/27/09	Brandi Stuart	Added Updates to Boiler Table
8/7/09	Brandi Stuart	Updated HFZ Table
1/11/10	Diane Garcia	Liability/Warranty sentence added at bottom of 1 <sup>st</sup> page
10/4/11	M. McCabe Janicki	Removed IFU sheet and placed in Forms
9/25/12	Edwena Jones	Added digital signature format
12/12/12	M. McCabe Janicki	Added an extra line between letter signature block and the word "Enclosure". Also, added a missing digit in 4-digit extension on letterhead zip code: "002" should be "0002".



# LUMENIS

## INDICATIONS FOR USE STATEMENT

510(k) Number (if known): K130164

Device Name: *FiberLase Endure C02 Fiber*

Indications for Use: The *FiberLase Endure C02 Fiber* is intended for use in surgical procedures requiring ablation, vaporization, excision, incision and coagulation of soft tissue. *FiberLase Endure C02 Fiber* is indicated for use in open surgical procedures such as ENT surgery and laparoscopy and endoscopic procedures.

Prescription Use:  (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use:  (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE -CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Neil R Ogden  
2013.03.22 16:12:33-04'00'

(Division Sign-Off) for MXM  
Division of Surgical Devices  
510(k) Number K130164

1 of 1

*FiberLase ENDURE C02 Fiber- 5/0k Notification*



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
Center for Devices and Radiological Health  
Document Control Center WO66-G609  
10903 New Hampshire Avenue  
Silver Spring, MD 20993-0002

March 11, 2013.

LUMENIS LTD.  
C/O QSITE  
31 HAAVODA ST.  
BINYAMINA  
ISRAEL 30500  
ATTN: YORAM LEVY

510k Number: K130164

Product: FIBERLASE ENDURE CO2 FIBER

The additional information you have submitted has been received.

We will notify you when the processing of this submission has been completed or if any additional information is required. Please remember that all correspondence concerning your submission **MUST** be sent to the Document Mail Center at the above letterhead address. Correspondence sent to any address other than the one above will not be considered as part of your official premarket notification submission. Also, please note the new Blue Book Memorandum regarding Fax and E-mail Policy entitled, "Fax and E-Mail Communication with Industry about Premarket Files Under Review. Please refer to this guidance for information on current fax and e-mail practices at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm089402.htm>. On August 12, 2005 CDRH issued the Guidance for Industry and FDA Staff: Format for Traditional and Abbreviated 510(k)s. This guidance can be found at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm084365.htm>. Please refer to this guidance for assistance on how to format an original submission for a Traditional or Abbreviated 510(k).

The Safe Medical Devices Act of 1990, signed on November 28, states that you may not place this device into commercial distribution until you receive a letter from FDA allowing you to do so. As in the past, we intend to complete our review as quickly as possible. Generally we do so in 90 days. However, the complexity of a submission or a requirement for additional information may occasionally cause the review to extend beyond 90 days. Thus, if you have not received a written decision or been contacted within 90 days of our receipt date you may want to check with FDA to determine the status of your submission.

**Please ensure that whether you submit a 510(k) Summary as per 21 CFR 807.92, or a 510(k) Statement as per 21 CFR 807.93, it meets the content and format regulatory requirements.**

If you have procedural questions, please contact the Division of Small Manufacturers International and Consumer Assistance (DSMICA) at (301)796-7100 or at their toll-free number (800)638-2041, or contact the 510k staff at (301)796-5640.

Sincerely,

510(k) Staff

**Pugh, Dominique \***

---

**From:** Pugh, Dominique \*  
**Sent:** Monday, March 11, 2013 6:30 PM  
**To:** Yoram@qsitemed.com  
**Cc:** DCCLetters  
**Subject:** K130164/S001 ACK LETTER  
**Attachments:** CrystalViewerCACMYO0K.rtf

**Tracking:**

**Recipient**

**Delivery**

Yoram@qsitemed.com

DCCLetters

Delivered: 3/11/2013 6:30 PM



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
Center for Devices and Radiological Health  
Document Control Center WO66-G609  
10903 New Hampshire Avenue  
Silver Spring, MD 20993-0002

January 30, 2013

LUMENIS LTD.  
C/O QSITE  
31 HAAVODA ST.  
BINYAMINA  
ISRAEL 30500  
ATTN: YORAM LEVY

510k Number: K130164

Received: 1/29/2013

Product: FIBERLASE ENDURE CO2 FIBER

The Food and Drug Administration (FDA), Center for Devices and Radiological Health (CDRH), has received the Premarket Notification, (510(k)), you submitted in accordance with Section 510(k) of the Federal Food, Drug, and Cosmetic Act(Act) for the above referenced product and for the above referenced 510(k) submitter. Please note, if the 510(k) submitter is incorrect, please notify the 510(k) Staff immediately. We have assigned your submission a unique 510(k) number that is cited above. Please refer prominently to this 510(k) number in all future correspondence that relates to this submission. We will notify you when the processing of your 510(k) has been completed or if any additional information is required. **YOU MAY NOT PLACE THIS DEVICE INTO COMMERCIAL DISTRIBUTION UNTIL YOU RECEIVE A LETTER FROM FDA ALLOWING YOU TO DO SO.**

Please remember that all correspondence concerning your submission **MUST** be sent to the Document Mail Center (DMC) at the above letterhead address. Correspondence sent to any address other than the one above will not be considered as part of your official 510(k) submission.

On September 27, 2007, the President signed an act reauthorizing medical device user fees for fiscal years 2008 - 2012. The legislation - the Medical Device User Fee Amendments of 2007 is part of a larger bill, the Food and Drug Amendments Act of 2007. Please visit our website at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/Overview/MedicalDeviceUserFeeandModernizationActMDUFMA/default.htm>

for more information regarding fees and FDA review goals. In addition, effective January 2, 2008, any firm that chooses to use a standard in the review of ANY new 510(k) needs to fill out the new standards form (Form 3654) and submit it with their 510(k). The form may be found at <http://www.fda.gov/AboutFDA/ReportsManualsForms/Forms/default.htm>.

We remind you that Title VIII of the Food and Drug Administration Amendments Act of 2007 (FDAAA) amended the PHS Act by adding new section 402(j) (42 U.S.C. § 282(j)), which expanded the current database known as ClinicalTrials.gov to include mandatory registration and reporting of results for applicable clinical trials of human drugs (including biological products) and devices. Section 402(j) requires that a certification form <http://www.fda.gov/AboutFDA/ReportsManualsForms/Forms/default.htm> accompany 510(k)/HDE/PMA submissions. The agency has issued a draft guidance titled: "Certifications To Accompany Drug, Biological

Product, and Device Applications/Submissions: Compliance with Section 402(j) of The Public Health Service Act, Added By Title VIII of The Food and Drug Administration Amendments Act of 2007”

<http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/HowtoMarketYourDevice/PremarketSubmissions/PremarketNotification510k/ucm134034.htm>. According to the draft guidance, 510(k) submissions that do not contain clinical data do not need the certification form.

Please note the following documents as they relate to 510(k) review: 1) Guidance for Industry and FDA Staff entitled, “Interactive Review for Medical Device Submissions: 510(k)s, Original PMAs, PMA Supplements, Original BLAs and BLA Supplements”. This guidance can be found at

<http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm089402.htm>.

Please refer to this guidance for information on a formalized interactive review process. 2) Guidance for Industry and FDA Staff entitled, "Format for Traditional and Abbreviated 510(k)s". This guidance can be found at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm084365.htm>. Please refer to this guidance for assistance on how to format an original submission for a Traditional or Abbreviated 510(k).

In all future premarket submissions, we encourage you to provide an electronic copy of your submission. By doing so, you will save FDA resources and may help reviewers navigate through longer documents more easily. Under CDRH's e-Copy Program, you may replace one paper copy of any premarket submission (e.g., 510(k), IDE, PMA, HDE) with an electronic copy. For more information about the program, including the formatting requirements, please visit our web site at

<http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/HowtoMarketYourDevice/PremarketSubmissions/ucm134508.html>. In addition, the 510(k) Program Video is now available for viewing on line at

<http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/HowtoMarketYourDevice/PremarketSubmissions/PremarketNotification510k/ucm070201.htm>.

**Please ensure that whether you submit a 510(k) Summary as per 21 CFR 807.92, or a 510(k) Statement as per 21 CFR 807.93, it meets the content and format regulatory requirements.**

Lastly, you should be familiar with the regulatory requirements for medical devices available at Device Advice <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/default.htm>. If you have questions on the status of your submission, please contact DSMICA at (301)796-7100 or the toll-free number (800)638-2041, or at their internet address <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/default.htm>. If you have procedural questions, please contact the 510(k) Staff at (301)796-5640.

Sincerely,

510(k) Staff

January 21, 2013  
Food and Drug Administration  
Center for Medical Devices and Radiological Health  
Office of Device Evaluation  
10903 New Hampshire Avenue  
White Oak, Bldg 66, Room 2534  
Silver Spring MD 20990  
USA

**Re: Traditional Premarket Notification for the Lumenis *FiberLase Endure* eCopy Cover Letter**

Dear Madam/Sir:

In accordance with the eCopy Program for Medical Device Submissions Lumenis, ("the Company") is submitting the enclosed Traditional 510(k) Notification (21 CFR 807.90(e)) eCopy that escorts the printed hard copy of the same submission at least ninety days before the Company intends to market in the United States the **FiberLase Endure** System.

**The eCopy is an exact duplicate of the paper copy.**

This letter is written following an ecopy hold letter from January 8, 2013.

As required by MDUFMA of 2002, a copy of the Medical Device User Fee Cover Sheet (Payment Identification Number **MD6066445-956733** is provided in Section 1 of this application.

Presented below is a summary of applicable information related to this submission. Detailed information is provided in Sections 1 to 16 of the submission.

**Applicant Name:**

Lumenis Ltd.  
13 Hayetzira Street  
Yokneam Industrial Park  
Yokneam 20692 Israel  
Tel. (972)4-959 9000  
Fax: (9724-959 9050

**Contact Person:**

Yoram Levy  
Haavoda 31,  
Binyamina 30500  
Israel  
Phone: (972)4-638-8837  
Fax: (972)4-638-0510  
E-mail: [Yoram@qsitemed.com](mailto:Yoram@qsitemed.com)

***Lumenis ResurFX 510k Notification Cover Letter***



**Trade Name:** Lumenis **FiberLase Endure**

Please direct any correspondence concerning this submission to Yoram Levy of QSite, Lumenis regulatory consultant, at [Yoram@qsite.com](mailto:Yoram@qsite.com), Tel (972)4-638-8837 or Fax (972)4-638-0510. Upon a finding of substantial equivalence, please send me a copy of the signed substantial equivalence letter by facsimile to (972)4-638-0510 or by Email.

Sincerely,

A handwritten signature in black ink that reads 'Yoram Levy'.

Lumenis Ltd. QA/RA Consultant

***Lumenis ResurFX 510k Notification Cover Letter***



## **Special 510(k) Submission**

# **FiberLase Endure CO2 Fiber**

Lumenis Ltd.  
13 Hayetzira Street  
Yokneam Industrial Park  
Yokneam 20692 Israel  
Tel. (972)4-959 9000  
Fax: (9724-959 9050





## **Special 510(k) PREMARKET NOTIFICATION** **CHECKLIST**

<b>Item</b>	<b>Comments</b>	
1.	510(k) Checklist	
2.	User fee Cover sheet Form	See 510(k) notice section A
3.	CDRH Submission Cover Sheet	See 510(k) notice section B
4.	Cover letter	See 510(k) notice section C
5.	PMA (Premarket Approval) cover page	NA – this is not a class III device
6.	PMN (Premarket Notification) checklist	See 510(k) notice page 3
7.	Table of contents	See 510(k) notice page i
8.	<b>Indication for Use Statement</b>	See 510(k) notice section 1
9.	<b>510(k) Summary</b>	See 510(k) notice section 2
10.	a) Administrative information (Applicant's Name, Contact Person ...)	See 510(k) notice section 2
11.	b) Name of the device (classification, Trade Name)	See 510(k) notice section 2
12.	c) Equivalence devices	See 510(k) notice section 2
13.	d) Description of the device	See 510(k) notice section 2
14.	e) Intended use	See 510(k) notice section 2
15.	f) Summary of the technological characteristics	See 510(k) notice section 2
16.	g) Performance data – non clinical	See 510(k) notice section 2
17.	h) Performance data – clinical	See 510(k) notice section 2
18.	i) Conclusions drawn from the non clinical and clinical	See 510(k) notice section 2
19.	<b>Truthful and Accuracy Statement</b>	See 510(k) notice section 3
20.	<b>Class III Summary and Certification</b>	NA - this is not a class III device
21.	<b>Financial Certification or Disclosure Statement</b>	See 510(k) notice section 4
22.	<b>Declaration of Conformity and Summary Reports</b>	See 510(k) notice section 5
23.	<b>Executive Summary</b>	See 510(k) notice section 6
24.	a) Description of the device, including the indications for use and technology	See 510(k) notice section 6
25.	b) Device comparison	See 510(k) notice section 6
26.	c) Summary for any performance testing in the	See 510(k) notice section 6



	submission	
27.	<b>Device description</b>	See 510(k) notice section 7
28.	a) General description	See 510(k) notice section 7
29.	b) Device performance specification	See 510(k) notice section 7
30.	c) Functional description	See 510(k) notice section 7
31.	d) HW description	See 510(k) notice section 7
32.	e) Device design requirements	See 510(k) notice section 7
33.	f) Diagrams, dimensions, schematics drawing and photos	See 510(k) notice section 7
34.	g) Patient contacting components and their respective materials	See 510(k) notice section 7
35.	<b>Substantial equivalence discussion</b>	See 510(k) notice section 8
36.	a) Predicate device, its trade name, model number, 510(k) submitter/holder, and 510(k) number	See 510(k) notice section 8
37.	b) Comparison discussion	See 510(k) notice section 8
38.	<b>Proposed Labeling</b>	See 510(k) notice section 9
39.	a) Device labels table	See 510(k) notice section 9
40.	b) Instructions Manual	See 510(k) notice section 9
41.	<b>Sterilization and Shelf Life</b>	See 510(k) notice section 10
42.	<b>Biocompatibility</b>	See 510(k) notice section 11
43.	a) Identification of material and manufacture as predicate device	See 510(k) notice section 11
44.	<b>Software</b>	See 510(k) notice section 12
45.	<b>Electromagnetic Compatibility and Electrical Safety</b>	See 510(k) notice section 13
46.	<b>Performance Testing - Bench</b>	See 510(k) notice section 14
47.	<b>Performance Testing - Animal</b>	See 510(k) notice section 15
48.	<b>Performance Testing - Clinical</b>	See 510(k) notice section 16
49.	<b>Other</b>	None

# TABLE OF CONTENTS

**SECTION A – MEDICAL DEVICE USER FEE COVER SHEET ..... 1**

**SECTION B – PREMARKET REVIEW SUBMISSION COVER SHEET ..... 1**

**SECTION C – 510K COVER LETTER ..... 1**

**SECTION 1 – INDICATION FOR USE STATEMENT ..... 1**

**SECTION 2- 510(K) SUMMARY..... 1**

**SECTION 3 – TRUTHFUL AND ACCURACY STATEMENT ..... 1**

**SECTION 4 – FINANCIAL CERTIFICATION OR DISCLOSURE STATEMENT ..... 1**

**SECTION 5 – DECLARATIONS OF CONFORMITY AND SUMMARY REPORTS..... 1**

**SECTION 6 –EXECUTIVE SUMMARY..... 1**

**SECTION 7 – DEVICE DESCRIPTION..... 1**

**7.1. GENERAL DESCRIPTION..... ERROR! BOOKMARK NOT DEFINED.**

**7.2. SYSTEM MAIN PARTS..... ERROR! BOOKMARK NOT DEFINED.**

**7.3. TECHNICAL PARAMETERS ..... ERROR! BOOKMARK NOT DEFINED.**

**SECTION 8 – SUBSTANTIAL EQUIVALENCE DISCUSSION ..... 1**

**8.1. PREDICATE DEVICES..... 2**

**8.2. COMPARISON DISCUSSION ..... 2**

**SECTION 9 - PROPOSED LABELING..... 1**

**9.1 INTENDED USE ..... 2**

**9.2 INSTRUCTIONS FOR USE ..... 2**

**9.3 DEVICE LABELS..... 2**

**9.4 PACKAGING LABELS ..... ERROR! BOOKMARK NOT DEFINED.**

**9.5 PACKAGE INSERT ..... ERROR! BOOKMARK NOT DEFINED.**

**9.6 QUICK GUIDES ..... ERROR! BOOKMARK NOT DEFINED.**

**SECTION 10- STERILIZATION AND SHELF LIFE..... 1**

**SECTION 11- BIOCOMPATIBILITY ..... 7**

**SECTION 12: SOFTWARE..... 2**

**SECTION 13: ELECTROMAGNETIC COMPATIBILITY AND ELECTRICAL SAFETY ..... 4**

**13.1 EMC ..... ERROR! BOOKMARK NOT DEFINED.**

**13.2 ELECTRICAL SAFETY ..... ERROR! BOOKMARK NOT DEFINED.**

**SECTION 14 - PERFORMANCE TESTING BENCH ..... 1**

**SECTION 15 - PERFORMANCE TESTING ANIMAL ..... 1**

**SECTION 16 - PERFORMANCE TESTING CLINICAL..... 1**

**ATTACHMENTS 1-13**

# **Section A – Medical Device User Fee Cover Sheet**



DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION <b>MEDICAL DEVICE USER FEE COVER SHEET</b>		PAYMENT IDENTIFICATION NUMBER: (b) (4) Write the Payment Identification number on your check.
A completed cover sheet must accompany each original application or supplement subject to fees. If payment is sent by U.S. mail or courier, please include a copy of this completed form with payment. Payment and mailing instructions can be found at: <a href="http://www.fda.gov/oc/mdufma/coversheet.html">http://www.fda.gov/oc/mdufma/coversheet.html</a>		
1. COMPANY NAME AND ADDRESS (include name, street address, city state, country, and post office code)  LUMENIS LTD New Ind. Park Yokneam 20692 IL  1.1 EMPLOYER IDENTIFICATION NUMBER (EIN)	2. CONTACT NAME Jace McLane  2.1 E-MAIL ADDRESS regulatory.affair@lumenis.com  2.2 TELEPHONE NUMBER (include Area code) 801-6562328  2.3 FACSIMILE (FAX) NUMBER (Include Area code)	
3. TYPE OF PREMARKET APPLICATION (Select one of the following in each column; if you are unsure, please refer to the application descriptions at the following web site: <a href="http://www.fda.gov/oc/mdufma">http://www.fda.gov/oc/mdufma</a> ) <u>Select an application type:</u> <input checked="" type="checkbox"/> Premarket notification(510(k)); except for third party <input type="checkbox"/> 513(g) Request for Information <input type="checkbox"/> Biologics License Application (BLA) <input type="checkbox"/> Premarket Approval Application (PMA) <input type="checkbox"/> Modular PMA <input type="checkbox"/> Product Development Protocol (PDP) <input type="checkbox"/> Premarket Report (PMR) <input type="checkbox"/> Annual Fee for Periodic Reporting (APR) <input type="checkbox"/> 30-Day Notice		
3.1 Select a center <input checked="" type="checkbox"/> CDRH <input type="checkbox"/> CBER 3.2 Select one of the types below <input checked="" type="checkbox"/> Original Application <u>Supplement Types:</u> <input type="checkbox"/> Efficacy (BLA) <input type="checkbox"/> Panel Track (PMA, PMR, PDP) <input type="checkbox"/> Real-Time (PMA, PMR, PDP) <input type="checkbox"/> 180-day (PMA, PMR, PDP)		
4. ARE YOU A SMALL BUSINESS? (See the instructions for more information on determining this status) <input type="checkbox"/> YES, I meet the small business criteria and have submitted the required qualifying documents to FDA <input checked="" type="checkbox"/> NO, I am not a small business  4.1 If Yes, please enter your Small Business Decision Number:		
5. FDA WILL NOT ACCEPT YOUR SUBMISSION IF YOUR COMPANY HAS NOT PAID AN ESTABLISHMENT REGISTRATION FEE THAT IS DUE TO FDA. HAS YOUR COMPANY PAID ALL ESTABLISHMENT REGISTRATION FEES THAT ARE DUE TO FDA? <input checked="" type="checkbox"/> YES (All of our establishments have registered and paid the fee, or this is our first device, and we will register and pay the fee within 30 days of FDA's approval/clearance of this device.) <input type="checkbox"/> NO (If "NO," FDA will not accept your submission until you have paid all fees due to FDA. This submission will not be processed; see <a href="http://www.fda.gov/cdrh/mdufma">http://www.fda.gov/cdrh/mdufma</a> for additional information)		
6. IS THIS PREMARKET APPLICATION COVERED BY ANY OF THE FOLLOWING USER FEE EXCEPTIONS? IF SO, CHECK THE APPLICABLE EXCEPTION. <input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates <input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population <input type="checkbox"/> This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only <input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially		
7. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION FOR WHICH FEES WERE WAIVED DUE TO SOLE USE IN A PEDIATRIC POPULATION THAT NOW PROPOSES CONDITION OF USE FOR ANY ADULT POPULATION? (If so, the application is subject to the fee that applies for an original premarket approval application (PMA)). <input type="checkbox"/> YES <input checked="" type="checkbox"/> NO		
PAPERWORK REDUCTION ACT STATEMENT Public reporting burden for this collection of information is estimated to average 18 minutes per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to the address below.  Department of Health and Human Services, Food and Drug Administration, Office of Chief Information Officer, 1350 Piccard Drive, 4th Floor Rockville, MD 20850 [Please do NOT return this form to the above address, except as it pertains to comments on the burden estimate.]		
8. USER FEE PAYMENT AMOUNT SUBMITTED FOR THIS PREMARKET APPLICATION (b) (4) <span style="float: right;">09-Jan-2013</span>		

Form FDA 3601 (01/2007)



# **Section B – Premarket Review Submission Cover Sheet**

B-1

*FiberLase ENDURE – 510k Notification*



## **Section C – 510k Cover Letter**



January 20, 2013

Food and Drug Administration

Center for Medical Devices and Radiological Health

Office of Device Evaluation

Document Mail Center (HFZ-401)

9200 Corporate Boulevard

Rockville, Maryland 20850 USA

**Re: Special Premarket Notification 510(k) Notification  
For *FiberLase Endure CO2 Fiber***

Dear Madam/Sir:

In accordance with section 510(k) of the Federal Food, Drug, and Cosmetic Act, and with the Code of Federal Regulations, 21 CFR 807, Subpart E, Lumenis Ltd. ("Lumenis" or "the Company"), hereby submits a special 510(k) Premarket Notification for *FiberLase Endure CO2 Fiber*.

For the convenience of reviewers, the Agency's "Screening Checklist for All Premarket Notifications" is included in this submission. The checklist identifies the sections within this submission where the required information is located.

As required by MDUFMA of 2002, a copy of the Medical Device User Fee Cover Sheet (Payment Identification Number (b) (4) [REDACTED]) is provided in Section A of this application.

Presented below is a summary of applicable information related to this submission. Detailed information is provided in Sections 1 to 16 of the submission.





**Applicant Name:** Lumenis Ltd.  
13 Hayetzira Street  
Yokneam Industrial Park  
Yokneam 20692 Israel  
Tel. (972)4-959-9000  
Fax: (972)4-959-9050

**Contact Person:** Yoram Levy, Qsite  
31 Haavoda St.  
Binyamina, ISRAEL  
Tel: (972)4-638-8837  
Fax: (972)4-638-0510

**Device Trade Name:** *FiberLase Endure CO2 Fiber*

**Device Type:** CO2 Laser fiber

**Classification:** **Regulation Name:** Laser surgical instrument for use in general and plastic surgery and in dermatology  
**Product Code:** GEX  
**Regulation No:** 21 CFR 878.4810  
**Class:** II  
**Classification Panel:** General & Plastic Surgery

**Indications for Use:** The *FiberLase Endure CO2 Fiber* is intended for use in surgical procedures requiring ablation, vaporization, excision, incision and coagulation of soft tissue. The *FiberLase Endure CO2 Fiber* is indicated for use in open surgical procedures such as ENT surgery and laparoscopy and endoscopic procedures.



**Predicate device:** **FiberLase CO2 Laser WaveGuide (K100384)**  
**Regulation Number:** 21 CFR 878.4810  
**Regulation Name:** Laser surgical instrument for use in  
general and plastic surgery and in dermatology  
**Regulatory Class:** II  
**Product Code:** GEX

**Type of submission:** Special 510k



**Design and Use of the Device:**

Question	YES	NO
Is the device intended for prescription use (21 CFR 801 Subpart D)? <sup>A</sup>	X	
Is the device intended for over-the-counter use (21 CFR 807 Subpart C)?		X
Does the device contain components derived from a tissue or other biologic source?		X
Is the device provided sterile?	X	
Is the device intended for single use?		X
Is the device a reprocessed single use device?		X
If yes, does this device type require reprocessed validation data?		N/A
Does the device contain a drug?		X
Does the device contain a biologic?		X
Does the device use software?		X
Does the submission include clinical information?		X
Is the device implanted?		X

**Confidentiality:**

(b) (4)

[Redacted]

[Redacted]

[Redacted]

C-5

*FiberLase Endure CO2 Fiber – 510k Notification*



(b) (4)

[Redacted content]

We trust that the information contained in this 510(k) Notification will be sufficient to enable FDA to find that the *FiberLase Endure CO2 Fiber* is substantially equivalent to its predicate device.

Please direct any correspondence concerning this submission to Yoram Levy of QSite, Lumenis regulatory consultant, at [Yoram@qsite.com](mailto:Yoram@qsite.com), Tel (972)4-638-8837 or Fax (972)4-638-0510. Upon a finding of substantial equivalence, please send me a copy of the signed substantial equivalence letter by facsimile to (972)4-638-0510 or by Email.

Sincerely,

  
Yoram Levy  
Lumenis Ltd. QA/RA Consultant

CC: Assaf Gelstein, Lumenis Ltd. Platform Manager



# **Section 1 – Indication for Use Statement**



**INDICATIONS FOR USE STATEMENT**

**510(k) Number (if known):**

**Device Name:** *FiberLase Endure CO2 Fiber*

**Indications for Use:** The *FiberLase Endure CO2 Fiber* is intended for use in surgical procedures requiring ablation, vaporization, excision, incision and coagulation of soft tissue. *FiberLase Endure CO2 Fiber* is indicated for use in open surgical procedures such as ENT surgery and laparoscopy and endoscopic procedures.

Prescription Use   X   AND/OR Over-The-Counter Use             
(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE -CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-off)  
Division of Anesthesiology, General Hospital, Infection Control and Dental Devices  
510(k) Number





## **Section 2- 510(k) Summary**





**510(K) SUMMARY**

***FiberLase Endure CO2 Fiber***

**510(k) Number K\_\_\_\_\_**

**Applicant's Name:** Lumenis Ltd.

13 Hayetzira Street  
Yokneam Industrial Park  
Yokneam 20692 Israel  
Tel. (972)4-959-9000  
Fax: (972)4-959-9050

**Contact Person:** Yoram Levy, Qsite

31 Haavoda St.  
Binyamina, Israel 30500  
Tel (972)4-638-8837  
Fax (972)4-638-0510  
Yoram@qsite.com

**Trade Name:** *FiberLase Endure CO2 Fiber*

**Device Type:** *CO2 Laser fiber*

**Preparation Date:** January 20, 2013

**Classification:** **Regulatory Name:** Laser surgical instrument for use in general and plastic surgery and in dermatology

**Product Code:** GEX

**Regulation No:** 21 CFR 878.4810

**Class:** II

**Classification Panel:** General & Plastic Surgery

**Device Description:**



The *FiberLase Endure CO2 Fiber* is a hollow, semi rigid, light-conducting delivery fiber, 2 meter length, designed to transmit laser energy from the CO2 laser system to the treatment site. The fiber also transmits a low power red diode or helium neon laser aiming beam to assist in targeting the tissue to be treated.

**Intended Use Statement:**

The *FiberLase Endure CO2 Fiber* is intended for use in surgical procedures requiring ablation, vaporization, excision, incision and coagulation of soft tissue. The *FiberLase Endure CO2 Fiber* is indicated for use in open surgical procedures such as ENT surgery and laparoscopy and endoscopic procedures.

**Predicate Devices:** Substantial equivalence to the following predicate device is claimed:

Device Name	510k No	Date of Clearance
FiberLase CO2 Laser WaveGuide	K100384	April 12 , 2010

**Performance Standards**

*FiberLase Endure CO2 Fiber* was tested and complies with the following standards:

- ANSI/AAMI/ISO 11135-1:2007 Sterilization of health care products — Ethylene oxide
- AAMI TIR30:2001 Product adoption and process equivalency for ethylene oxide sterilization
- ISO 14971-1:2007 Risk management for medical devices
- IEC 60601-2-22 ed3.0:2007 – Medical Electrical Equipment - Part 2-22: Particular Requirements for Basic Safety and Essential Performance of Surgical, Cosmetic, Therapeutic and Diagnostic Laser Equipment.
- IEC 60825-1:2007 – Safety of Laser Products - Part 1: Equipment Classification, Requirements and User's Guide
- ANSI/AAMI/ISO 17665-1:2006 Sterilization of health care products — Moist Heat.



A detailed description follows in **Section 14**.

### **Performance Testing**

Performance testing demonstrated that the *FiberLase Endure CO2 Fiber* is as safe and effective as the cleared predicate device.

### **Comparison with the Predicate Devices**

The *FiberLase Endure CO2 Fiber* is a modification to its predicate device, the FDA-cleared Lumenis Fiberlase CO2 laser WaveGuide (K100384).

The intended use of the *FiberLase Endure CO2 Fiber* is identical to the intended use of its predicate.

Both the *FiberLase Endure CO2 Fiber* and the Lumenis Fiberlase CO2 laser WaveGuide systems are fibers that transmit laser energy from the laser system to the treatment site. Both devices are comprised of a laser connector and a 2 meter long delivery fiber.

The structures, the materials and the dimensions of the *FiberLase Endure CO2 Fiber* are identical to the cleared Lumenis Fiberlase CO2 laser WaveGuide fiber.

The minor difference between the *FiberLase Endure CO2 Fiber* System and its predicate device is an addition of an Autoclave sterilization method for reprocessing the fiber instead of single use ETO. The minor differences do not raise any new questions of safety or efficacy. Moreover, performance testing demonstrated that the *FiberLase Endure CO2 Fiber* is as safe and effective as the predicate device. Thus, the *FiberLase Endure CO2 Fiber* is substantially equivalent to Lumenis FiberLase CO2 Laser WaveGuide (K100384).



## **Section 3 – Truthful and Accuracy Statement**



**PREMARKET NOTIFICATION**

**TRUTHFUL AND ACCURATE STATEMENT\***

**(As Required by 21 CFR 807.87 (j))**

I certify that, in my capacity as Platform Manager at Lumenis, I believe to the best of my knowledge, that all data and information submitted in the premarket notification are truthful and accurate and that no material fact has been omitted.

A handwritten signature in black ink, appearing to read 'Assaf Gelstein', written over a horizontal line.

**Signature**

Assaf Gelstein

**Typed Name and Title**

Lumenis Ltd.

**Company**

January 20, 2013

**Date**

---

**Premarket Notification 510(k) Number**



## **Section 4 – Financial Certification or Disclosure Statement**



**SECTION 4– FINANCIAL CERTIFICATION OR DISCLOSURE**  
**STATEMENT**

No clinical study was conducted. Therefore, no Financial Certification or Disclosure Statement is needed.



# **Section 5 – Declarations of Conformity and Summary Reports**





## COMPLIANCE WITH STANDARDS

The *FiberLase Endure CO2 Fiber* complies with:

- IEC 60601-2-22 ed3.0:2007 – Medical Electrical Equipment - Part 2-22: Particular Requirements for Basic Safety and Essential Performance of Surgical, Cosmetic, Therapeutic and Diagnostic Laser Equipment.
- IEC 60825-1:2007 – Safety of Laser Products - Part 1: Equipment Classification, Requirements and User's Guide

When combines with the laser system

In addition, the *FiberLase Endure CO2 Fiber* complies with the following voluntary standards (certification for compliance with the following standards follows in this section):

1. ANSI/AAMI/ISO 11135-1:2007 Sterilization of health care products — Ethylene oxide.
2. ANSI/AAMI/ISO 17665-1:2006 Sterilization of health care products — Moist Heat.
3. AAMI TIR30:2011 A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices
4. AAMI TIR12:2010 Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers
5. ISO 14971-1:2007 Risk management for medical devices.



## VALIDATION AND VERIFICATION TESTING

A validation and verification testing report for the *FiberLase Endure CO2 Fiber* is attached to this submission. Testing result shows that the *FiberLase Endure CO2 Fiber* is safe and that it performs according to its specifications.



## STATEMENT OF CONFORMANCE TO PERFORMANCE STANDARDS

The company certifies that during the design process of *FiberLase Endure CO2 Fiber*, it was tested according to the performance standards.

The ETO sterilization validation was done according to ANSI/AAMI/ISO 11135-1:2007 Sterilization of health care products — Ethylene oxide and the Autoclave sterilization validation was done according to ANSI/AAMI/ISO 17665

The verification testing demonstrates that the functional requirement and system specifications were met.

The company also certifies that during production of *FiberLase Endure CO2 Fiber*, it will be validated for safety and integrity.

Signature:

A handwritten signature in black ink, appearing to read 'Ya'acov Yariv', with a long horizontal stroke extending to the right.

Date: January 20, 2013

Ya'acov Yariv

QA Manager

Lumenis



## DECLARATION OF CONFORMITY

(Traditional 510k)

To the best of my knowledge, the verification activities, as required by the risk analysis, for the *FiberLase Endure CO2 Fiber*, were performed by the designated individual(s) and the results demonstrated that the pre-determined acceptance criteria were met.

We, Lumenis, 13 Hayetzira Street 2006, The Industrial Park, Yokneam 20692 Israel, declare that our facility is in conformance with the design control requirements as specified in 21 CFR 820.30 and the records are available for review.

We further declare that qualified company employees performed all verification and validation activities that were required by the risk management file and the results demonstrated that the predetermined acceptance criteria were met.

A handwritten signature in black ink, appearing to read 'Yariv', written over a horizontal line.

Ya'acov Yariv  
QA Manager  
Lumenis

Date: January 20, 2013



## **Section 6 –Executive Summary**

6-1

*FiberLase Endure CO2 Fiber – 510k Notification*



## **SECTION 6 – EXECUTIVE SUMMARY**

### **6.1. *DEVICE DESCRIPTION***

The *FiberLase Endure CO2 Fiber*, is comprised of a laser connector and a 2 meter long delivery fiber. The fiber transmits laser energy from the laser system to the treatment site. The fiber also transmits a low-level red diode laser or helium-neon laser aiming beam to assist in targeting the tissue to be treated.

The *FiberLase Endure CO2 Fiber*, is a hollow, semi rigid, light-conducting tube designed to operate at the wavelengths of 10.6  $\mu\text{m}$  and 0.635  $\mu\text{m}$ . These wavelengths correspond to the output of the CO2 laser device and the red-diode or helium-neon aiming beam.

The *FiberLase Endure CO2 Fiber*, is a modification to the market-cleared Lumenis FiberLase CO2 Laser WaveGuide, which was cleared on April 12 2010, under K100384. The primary difference between the modified *FiberLase Endure CO2 Fiber*, and the market-cleared Lumenis FiberLase CO2 Laser WaveGuide is:

- Changing the fiber from single- to multi-use by Addition of an Autoclave sterilization method for reprocessing the fiber

The modified *FiberLase Endure CO2 Fiber* has undergone changes that do not affect the fiber's intended use. In addition, the changes were validated to demonstrate that the modified fiber is substantially equivalent to the cleared fiber and does not raise any new safety and the effectiveness questions.

### **6.2. *INDICATIONS FOR USE***

The *FiberLase Endure CO2 Fiber*, is intended for use in surgical procedures requiring ablation, vaporization, excision, incision and coagulation of soft tissue. The *FiberLase Endure CO2 Fiber*, is indicated for use in open surgical procedures such as ENT surgery and laparoscopy and endoscopic procedures.



### 6.3. PREDICATE DEVICE COMPARISON

The *FiberLase Endure CO2 Fiber*, is substantially equivalent to the following predicate device:

Device Name	510k No	Date of Clearance
Lumenis FiberLase CO2 Laser WaveGuide	K100384	April 12 2010

#### 6.3.1. Comparison Discussion:

The *FiberLase Endure CO2 Fiber*, is a modification to its predicate device, the FDA-cleared Lumenis Fiberlase CO2 laser WaveGuide (K100384).

The *FiberLase Endure CO2 Fiber* has same intended use as the cleared Lumenis FiberLase CO2 Laser WaveGuide (K100384).

The *FiberLase Endure CO2 Fiber*, System utilizes the same technology as the cleared Lumenis FiberLase CO2 Laser WaveGuide (K100384).

The fiber of the proposed system is made of the same material as the fiber of the cleared Lumenis FiberLase CO2 Laser WaveGuide (K100384).

The primary difference between the proposed *FiberLase Endure CO2 Fiber* and the market-cleared Lumenis FiberLase CO2 Laser WaveGuide is:

- Changing the fiber from single- to multi-use by Addition of an Autoclave sterilization method for reprocessing the fiber.

### 6.4. PERFORMANCE TESTING

(b) (4)



(b) (4)

A large rectangular area of the document is completely redacted with a solid grey fill.

## **6.5. CONCLUSION**

(b) (4)

A large rectangular area of the document is completely redacted with a solid grey fill.





**ATTACHMENTS TO SECTION 6**

---

**Attachment No. 1: Substantial Equivalence Comparison Table**



## **Section 7 – Device Description**



## **SECTION 7 – DEVICE DESCRIPTION**

### **7.1. GENERAL DESCRIPTION**

The *FiberLase Endure CO2 Fiber* is intended for use in surgical procedures requiring ablation, vaporization, excision, incision and coagulation of soft tissue. The *FiberLase Endure CO2 Fiber* is indicated for use in open surgical procedures such as ENT surgery and laparoscopy and endoscopic procedures.

The *FiberLase Endure CO2 Fiber* is comprised of a laser connector and a 2 meter long delivery reusable fiber. The fiber transmits laser energy from the laser system through the fiber to the treatment site. The fiber also transmits a low-level red diode or helium-neon aiming beam to assist in targeting the tissue to be treated.

The *FiberLase Endure CO2 Fiber* is a hollow, semi rigid, light-conducting tube designed to operate at the wavelengths of 10.6  $\mu\text{m}$  and 0.635  $\mu\text{m}$ . These wavelengths correspond to the output of the CO2 laser device and the red-diode or helium-neon aiming beam.

The *FiberLase Endure CO2 Fiber* is a modification to the market-cleared Lumenis FiberLase CO2 Laser WaveGuide, which was cleared on April 12 2010, under K100384. The primary difference between the modified **FiberLase Endure CO2 Fiber** and the market-cleared Lumenis FiberLase CO2 Laser WaveGuide is:

- Changing the fiber from single- to multi-use by Addition of an Autoclave sterilization method for reprocessing the fiber

The *FiberLase Endure CO2 Fiber* has undergone changes that do not affect the fiber's intended use. In addition, the changes were validated to demonstrate that the modified fiber is substantially equivalent to the original fiber and the changes do not affect the safety and the effectiveness of the fiber.



**7.2. FIBERLASE ENDURE CO2 FIBER**

(b) (4)

A large, solid grey rectangular redaction box covers the majority of the page content, starting below the section header and ending above the footer.



(b) (4)

A large, solid grey rectangular redaction box covering the majority of the page's content.

### **7.3. PRINCIPLE OF OPERATION**

(b) (4)

A large, solid grey rectangular redaction box covering the content of section 7.3.

### **7.4. FIBERLASE CO<sub>2</sub> LASER FIBERS-ACCESSORIES**

(b) (4)

A large, solid grey rectangular redaction box covering the content of section 7.4.



(b) (4)

A large, solid grey rectangular redaction box covers the majority of the page's content, starting below the LUMENIS logo and extending nearly to the bottom of the page.



**ATTACHMENTS TO SECTION 7**

---

**Attachment No. 2a: FiberLase Endure CO2 Fiber System User Manual**

**Attachment No. 2b: FiberLase Handpieces User Manual**

**Attachment No. 2c: FiberLase Fiber Renewal Kit User Manual**



# **Section 8 – Substantial Equivalence Discussion**





## SECTION 8 - SUBSTANTIAL EQUIVALENCE DISCUSSION

### 8.1. PREDICATE DEVICES

The *FiberLase Endure CO2 Fiber* is a reusable CO2 fiber, designed to transmit laser energy from a laser system through the fiber to the treatment site.

The *FiberLase Endure CO2 Fiber* is substantially equivalent to the following predicate device:

Device Name	510k No	Date of Clearance
Lumenis FiberLase CO2 Laser WaveGuide	K100384	April 12 2010

### 8.2. COMPARISON TABLE

Substantial equivalence table, which summarizes the similarities and differences between the *FiberLase Endure CO2 Fiber* and its predicate device, the Lumenis FiberLase CO2 Laser WaveGuide, is attached to this submission as Attachment no. 1 and follows in this section in table 1:

Characteristics	Proposed FiberLase Endure CO2 Fiber	Lumenis FiberLase CO2 Laser WaveGuide (K100384)
Intended use	(b) (4)	Same
Product Code		Same
Regulation No.		Same



Characteristics	Proposed FiberLase Endure CO2 Fiber	Lumenis FiberLase CO2 Laser WaveGuide (K100384)
<b>Classification</b>	(b) (4)	Same
<b>Operation</b>		Same
<b>Basic Structure</b>		Same
<b>Fiber dimensions</b>		Same. (other sizes are available)
<b>Materials</b>		Same
<b>Use</b>		Single Use
<b>Method of Sterilization</b>		The fiber is provided sterilized by Ethylene Oxide.
<b>Treatment's laser</b>		Same
<b>Aiming beam</b>		633 nm (Helium-Neon aiming beam)
<b>Maximum Input Power</b>		Same
<b>Performance Tests</b>		Non-clinical Laboratory performance testing was performed for the fiber functionality (e.g. Energy transmission)



Characteristics	Proposed FiberLase Endure CO2 Fiber	Lumenis FiberLase CO2 Laser WaveGuide (K100384)
Shelf Life	(b) (4)	Same
Packaging	(b) (4)	Same

**Table No. 1: Comparison between the Proposed *FiberLase Endure CO2 Fiber* and the FDA cleared Lumenis FiberLase CO2 Laser WaveGuide (K100384)**

**8.3. COMPARISON DISCUSSION**

(b) (4)





(b) (4)

A large, solid grey rectangular redaction box covering the majority of the page's content.

#### **8.4. CONCLUSION**

(b) (4)

A large, solid grey rectangular redaction box covering the majority of the page's content.



**ATTACHMENTS TO SECTION 8**

---

**Attachment No. 3: Lumenis FiberLase CO2 Laser WaveGuide (K100384)**

**510(k) Summary**



## **Section 9 - Proposed Labeling**



## **SECTION 9 – PROPOSED LABELING**

### ***9.1 INTENDED USE***

The *Lumenis FiberLase Endure CO2 Fiber* is intended for use in surgical procedures requiring ablation, vaporization, excision, incision and coagulation of soft tissue. The *FiberLase Endure CO2 Fiber* is indicated for use in open surgical procedures such as ENT surgery and laparoscopy and endoscopic procedures.

### ***9.2 INSTRUCTIONS FOR USE***

*FiberLase Endure CO2 Fiber's* User Manual is attached to this submission (Attachment No 2a). Also attached to this submission are *FiberLase Endure CO2 Fiber's* Handpieces User Manual (attachment No. 2b) and *FiberLase Endure CO2 Fiber's* Renewal Kit User Manual (attachment No. 2c).

### ***9.3 DEVICE LABELS***

The following describes the labeling of the *FiberLase Endure CO2 Fiber* (Figure No. 9-1):



**LUMENIS**<sup>®</sup>  
Enhancing Life. Advancing Technology.

**FiberLase™ ENDURE™**  
**CO<sub>2</sub> Fiber**

<b>REF</b>	AC-1148110		<b>Qty: 1</b>
<b>LOT</b>	4250312		2017-12

LENGTH 2 m    OD 1.04 mm

---

  
\*+AC-1148110%\*

  
\*+\$LUM112174250113%.\*

---

**STERILE** **EO**      

 **CONSULT INSTRUCTIONS FOR USE**       **DO NOT USE IF PACKAGE IS DAMAGED**

**Rx Only**      **CAUTION:**  
Federal law restricts this device to sale by  
or on the order of a physician

 **LUMENIS LTD.**  
P.O. Box 240  
Yokneam 20692, ISRAEL  
Tel.: +972 (4) 959-8000

 **Authorized Representative:**  
Lumenis (Germany) GmbH  
Heinrich-Hertz-Str. 3  
D-43303 Dreieich-Dreieichenhain  
Germany  
Tel: +49 (0)6103 - 8335-0  
Fax: +49 (0)6103 - 8335-300

**US Customer Support:**  
Lumenis Inc., USA  
Tel: +1 (877) 586-3647

LB-10000040

**FiberLase™ ENDURE™ CO<sub>2</sub> Fiber**

**Surgery Cycle –**

1	2	3	4	5
---	---	---	---	---









Lot#: 55550910 LB-10000050

*Figure 9-1: FiberLase Endure CO2 Fiber Labeling*





The following describes the labeling symbols explanations of the *FiberLase Endure CO2 Fiber* (Figure No. 9-2):

	Follow instructions for use
	CE Compliance
	Part number
	Lot number
	Expiration date, in yyyy-mm format
	EO sterilized
	Manufacturer's details
	Authorized representative
<b>Rx ONLY</b>	CAUTION: FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN. U.S. sales restriction label

*Figure 9-2: FiberLase Endure CO2 symbol explanations*



**ATTACHMENTS TO SECTION 9**

---

**None**



## **Section 10- Sterilization and Shelf Life**

10-1

*FiberLase ENDURE CO2 Fiber – 510k Notification*

**SECTION 10 – STERILIZATION AND SHELF LIFE**

**10.1 STERILIZATION AND CLEANING**

(b) (4)





(b) (4)

A large, solid grey rectangular redaction box covering the majority of the page's content.

**10.1.2. Fiber accessories sterilization and cleaning**

(b) (4)

A large, solid grey rectangular redaction box covering the majority of the page's content.



## **10.2 PACKAGING**

The *FiberLase Endure CO2 Fiber* is packed in a double Tyvek pouch, suitable for the ETO sterilization process, exactly as its predicate, the market-cleared Lumenis FiberLase CO2 laser WaveGuide (K100384).

The pouched *FiberLase Endure CO2 Fiber* is sent to the customer in a cardboard box. An envelope with the operating manual and warranty card is also enclosed in the main box. The *FiberLase Endure CO2 Fiber* packaging has been validated as suitable for shipping in accordance with ASTM 4169-96 "standard practice for performance testing of shipping containers and system" and other standards.

After each use, the fiber and the accessories can be reused, following manual cleaning and sterilization process. The autoclave sterilization process should be performed in a pouch designed for steam sterilization or in a Lumenis' sterilization tray designed for steam sterilization (supply with the fiber), wrapped with a medical grade steam sterilization wrap (According to ISO 17665-1).

## **10.3. SHELF LIFE**

(b) (4)

A large rectangular area of the document is completely redacted with a solid grey fill, obscuring all text and graphics that would otherwise be present in this section.



(b) (4)

A large, solid grey rectangular redaction box covers the majority of the page's content, starting below the header and ending above the footer.

**ATTACHMENTS TO SECTION 10**

---

**Attachment No. 4a: Reusable fiber cleaning validation report**

**Attachment No. 4b: Reusable fiber Autoclave sterilization validation report**

**Attachment No. 4c: Reusable FiberLase Handpieces and associated devices cleaning validation report**

**Attachment No. 4d: Reusable FiberLase Handpieces and associated devices Autoclave sterilization validation report**





# **Section 11- Biocompatibility**



**SECTION 11 – BIOCOMPTABILITY**  
**MATERIALS AND BIOCOMPATIBILITY**

The *FiberLase Endure CO2 Fiber* is considered as non-contact item as its predicate market-cleared Lumenis FiberLase CO2 laser WaveGuide (K100384).

In addition, the *FiberLase Endure CO2 Fiber* is made of the exactly the same materials as its predicate:

(b) (4)





**ATTACHMENTS TO SECTION 11**

---

**None**



## Section 12- Software



## **SECTION 12: SOFTWARE**

The *FiberLase Endure CO2 Fiber System* does not include software.



# **Section 13 - Electromagnetic Compatibility and Electrical Safety**



**SECTION 13: ELECTROMAGNETIC COMPATIBILITY AND  
ELECTRICAL SAFETY**

The *FiberLase Endure CO2 Fiber System* does not include Electrical and Electromagnetic systems.



**ATTACHMENTS TO SECTION 13**

---

**None**





# Section 14 - Performance Testing Bench



## **SECTION 14 – PERFORMANCE TESTING - BENCH**

### **14.1. Risk analysis**

(b) (4)

A large rectangular area of the document is redacted with a solid grey fill, covering the entire content of section 14.1.

### **14.2. Performance test**

(b) (4)

A large rectangular area of the document is redacted with a solid grey fill, covering the entire content of section 14.2.



(b) (4)

A large, solid grey rectangular redaction box covers the majority of the page's content, starting below the LUMENIS logo and extending nearly to the bottom of the page.



## **ATTACHMENTS TO SECTION 14**

---

**Attachment No. 5: FiberLase Endure CO2 Fiber Risk Analysis**

**Attachment No. 6: FiberLase Endure CO2 Fiber Functionality test report**

**Attachments No. 7: FDA Forms 3654 (Standards)**



# **Section 15 - Performance Testing Animal**



## **SECTION 15 – PERFORMANCE TESTING - ANIMAL**

The safety and efficacy of the *FiberLase Endure CO2 Fiber* was shown by bench testing and performance studies.

Lumenis believes that animal study is not needed to claim safety and efficacy of the device.



## **ATTACHMENTS TO SECTION 15**

---

**None**



# **Section 16 - Performance Testing Clinical**





## **SECTION 16 – PERFORMANCE TESTING - CLINICAL**

The safety and efficacy of the *FiberLase Endure CO2 Fiber* was shown by bench testing and performance studies.

Lumenis believes that clinical study is not needed to claim safety and efficacy of the device.



## **ATTACHMENTS TO SECTION 16**

---

### **Attachment No. 8: FDA Form 3674 Requirements of ClinicalTrials.gov**



## **List of Attachments**

Attachment No. 1: Substantial Equivalence Comparison Table

Attachment No. 2a: FiberLase Endure CO2 Fiber System User Manual

Attachment No. 2b: FiberLase Handpieces User Manual

Attachment No. 2c: FiberLase Fiber Renewal Kit User Manual

Attachment No. 3: Lumenis FiberLase CO2 Laser WaveGuide (K100384) 510(k) Summary

Attachment No. 4a: Cleaning validation reusable CO2 FiberLase Fibers

Attachment No. 4b: Sterilization validation for reusable CO2 FiberLase Fibers

Attachment No. 4c: Cleaning validation GLP report- Reusable FiberLase Handpieces and  
associated devices

Attachment No. 4d: Sterilization validation GLP report - Reusable FiberLase Handpieces  
and associated devices Autoclave - round 1

Attachment No. 4e: sterilization validation GLP report - Reusable FiberLase Handpieces  
and associated devices Autoclave - round 2

Attachment No. 5: FiberLase Endure CO2 Fiber Risk Analysis

Attachment No. 6: FiberLase Endure CO2 Fiber Functionality test report

Attachment No. 7: FDA Forms 3654 (Standards)

Attachment No. 8: FDA Form 3674 Requirements of ClinicalTrials.gov

Attachments



**Attachment No. 1:  
Substantial Equivalence Comparison Table**

Attachments



**Attachment No. 2a:  
FiberLase Endure CO2 Fiber System User Manual**

Attachments



**Attachment No. 2b:  
FiberLase Handpieces User Manual**

Attachments

*FiberLase Endure CO2 Fiber - 510k Notification*

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118



**Attachment No. 2c:  
FiberLase Fiber Renewal Kit User Manual**

Attachments

*FiberLase Endure CO2 Fiber - 510k Notification*

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118



**Attachment No. 3:  
Lumenis FiberLase CO2 Laser WaveGuide (K100384) 510(k)  
Summary**

Attachments

*FiberLase Endure CO2 Fiber - 510k Notification*

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOI STATUS@fda.hhs.gov or 301-796-8118





**Attachment No. 4a:  
Cleaning validation reusable CO2 FiberLase Fibers**

Attachments

*FiberLase Endure CO2 Fiber - 510k Notification*

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118



**Attachment No. 4b:  
Sterilization validation for reusable CO2 FiberLase Fibers**

Attachments

*FiberLase Endure CO2 Fiber - 510k Notification*

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOI STATUS@fda.hhs.gov or 301-796-8118



**Attachment No. 4c:  
Cleaning validation GLP report- Reusable FiberLase  
Handpieces and associated devices**

Attachments

*FiberLase Endure CO2 Fiber - 510k Notification*

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118



**Attachment No. 4d:  
Sterilization validation GLP report - Reusable FiberLase  
Handpieces and associated devices Autoclave - round I**

Attachments

*FiberLase Endure CO2 Fiber - 510k Notification*

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118



**Attachment No. 4e:  
Sterilization validation GLP report - Reusable FiberLase  
Handpieces and associated devices Autoclave - round 2**

Attachments

*FiberLase Endure CO2 Fiber - 510k Notification*

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118



**Attachment No. 5:  
FiberLase Endure CO2 Fiber Risk Analysis**

Attachments

*FiberLase Endure CO2 Fiber - 510k Notification*

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118



**Attachment No. 6:  
FiberLase Endure CO2 Fiber Functionality test report**

Attachments



**Attachment No. 7**  
**FDA Forms 3654 (Standards)**

Attachments

*FiberLase Endure CO2 Fiber - 510k Notification*

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118





## **Attachment No. 8**

### **FDA Form 3674 Requirements of ClinicalTrials.gov**

Attachments

*FiberLase Endure CO2 Fiber - 510k Notification*

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

DEPARTMENT OF HEALTH AND HUMAN SERVICES  
FOOD AND DRUG ADMINISTRATIONForm Approval  
OMB No. 9010-0120  
Expiration Date: May 31, 2007.  
See OMB Statement on page 5.**CDRH PREMARKET REVIEW SUBMISSION COVER SHEET**Date of Submission  
January 20, 2013User Fee Payment ID Number  
(b) (4)

FDA Submission Document Number (if known)

**SECTION A****TYPE OF SUBMISSION**

<b>PMA</b>	<b>PMA &amp; HDE Supplement</b>	<b>PDP</b>	<b>510(k)</b>	<b>Meeting</b>
<input type="checkbox"/> Original Submission <input type="checkbox"/> Premarket Report <input type="checkbox"/> Modular Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment <input type="checkbox"/> Licensing Agreement	<input type="checkbox"/> Regular (180 day) <input type="checkbox"/> Special <input type="checkbox"/> Panel Track (PMA Only) <input type="checkbox"/> 30-day Supplement <input type="checkbox"/> 30-day Notice <input type="checkbox"/> 135-day Supplement <input type="checkbox"/> Real-time Review <input type="checkbox"/> Amendment to PMA &HDE Supplement <input type="checkbox"/> Other	<input type="checkbox"/> Original PDP <input type="checkbox"/> Notice of Completion <input type="checkbox"/> Amendment to PDP	<input checked="" type="checkbox"/> Original Submission: <input type="checkbox"/> Traditional <input checked="" type="checkbox"/> Special <input type="checkbox"/> Abbreviated (Complete section I, Page 5) <input type="checkbox"/> Additional Information <input type="checkbox"/> Third Party	<input type="checkbox"/> Pre-510(K) Meeting <input type="checkbox"/> Pre-IDE Meeting <input type="checkbox"/> Pre-PMA Meeting <input type="checkbox"/> Pre-PDP Meeting <input type="checkbox"/> Day 100 Meeting <input type="checkbox"/> Agreement Meeting <input type="checkbox"/> Determination Meeting <input type="checkbox"/> Other (specify):
<b>IDE</b>	<b>Humanitarian Device Exemption (HDE)</b>	<b>Class II Exemption Petition</b>	<b>Evaluation of Automatic Class III Designation (De Novo)</b>	<b>Other Submission</b>
<input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement	<input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment	<input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information	<input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information	<input type="checkbox"/> 513(g) <input type="checkbox"/> Other (describe submission):

Have you used or cited Standards in your submission?  Yes  No (If Yes, please complete Section I, Page 5)**SECTION B****SUBMITTER, APPLICANT OR SPONSOR**

Company / Institution Name Lumenis Ltd.	Establishment Registration Number (if known) (b) (4)		
Division Name (if applicable)	Phone Number (including area code) (972)4-959 9000		
Street Address 13 Hayetzira Street	FAX Number (including area code) (9724-959 9050		
City Yokneam	State / Province	ZIP/Postal Code 20692	Country Israel
Contact Name Yoram Levy			
Contact Title Lumenis QA/RA Consultant		Contact E-mail Address Yoram@qsitemed.com	

**SECTION C****APPLICATION CORRESPONDENT (e.g., consultant, if different from above)**

Company / Institution Name Qsite	Establishment Registration Number (if known) (b) (4)		
Division Name (if applicable)	Phone Number (including area code) ( 972)4-638-8837		
Street Address 31 Haavoda St.	FAX Number (including area code) ( 972)4-638-0510		
City Binyamina	State / Province	ZIP/Postal Code 30500	Country Israel
Contact Name Yoram Levy			
Contact Title Qsite General Manager		Contact E-mail Address yoram@qsitemed.com	

SECTION D1			REASON FOR APPLICATION - PMA, PDP, OR HDE		
<input type="checkbox"/> Withdrawal <input type="checkbox"/> Additional or Expanded Indications <input type="checkbox"/> Request for Extension <input type="checkbox"/> Post-approval Study Protocol <input type="checkbox"/> Request for Applicant Hold <input type="checkbox"/> Request for Removal of Applicant Hold <input type="checkbox"/> Request to Remove or Add Manufacturing Site	<input type="checkbox"/> Change in design, component, or specification: <input type="checkbox"/> Software / Hardware <input type="checkbox"/> Color Additive <input type="checkbox"/> Material <input type="checkbox"/> Specifications <input type="checkbox"/> Other ( <i>specify below</i> )	<input type="checkbox"/> Location change: <input type="checkbox"/> Manufacturer <input type="checkbox"/> Sterilizer <input type="checkbox"/> Packager			
<input type="checkbox"/> Process change: <input type="checkbox"/> Manufacturing <input type="checkbox"/> Sterilization <input type="checkbox"/> Packaging <input type="checkbox"/> Other ( <i>specify below</i> )	<input type="checkbox"/> Labeling change: <input type="checkbox"/> Indications <input type="checkbox"/> Instructions <input type="checkbox"/> Performance <input type="checkbox"/> Shelf Life <input type="checkbox"/> Trade Name <input type="checkbox"/> Other ( <i>specify below</i> )	<input type="checkbox"/> Report Submission: <input type="checkbox"/> Annual or Periodic <input type="checkbox"/> Post-approval Study <input type="checkbox"/> Adverse Reaction <input type="checkbox"/> Device Defect <input type="checkbox"/> Amendment			
<input type="checkbox"/> Response to FDA correspondence:		<input type="checkbox"/> Change in Ownership <input type="checkbox"/> Change in Correspondent <input type="checkbox"/> Change of Applicant Address			
<input type="checkbox"/> Other Reason ( <i>specify</i> )					

SECTION D2			REASON FOR APPLICATION - IDE		
<input type="checkbox"/> New Device <input type="checkbox"/> New Indication <input type="checkbox"/> Addition of Institution <input type="checkbox"/> Expansion / Extension of Study <input type="checkbox"/> IRB Certification <input type="checkbox"/> Termination of Study <input type="checkbox"/> Withdrawal of Application <input type="checkbox"/> Unanticipated Adverse Effect <input type="checkbox"/> Notification of Emergency Use <input type="checkbox"/> Compassionate Use Request <input type="checkbox"/> Treatment IDE <input type="checkbox"/> Continued Access	<input type="checkbox"/> Change in: <input type="checkbox"/> Correspondent / Applicant <input type="checkbox"/> Design / Device <input type="checkbox"/> Informed Consent <input type="checkbox"/> Manufacturer <input type="checkbox"/> Manufacturing Process <input type="checkbox"/> Protocol - Feasibility <input type="checkbox"/> Protocol - Other <input type="checkbox"/> Sponsor	<input type="checkbox"/> Repose to FDA Letter Concerning: <input type="checkbox"/> Conditional Approval <input type="checkbox"/> Deemed Approved <input type="checkbox"/> Deficient Final Report <input type="checkbox"/> Deficient Progress Report <input type="checkbox"/> Deficient Investigator Report <input type="checkbox"/> Disapproval <input type="checkbox"/> Request Extension of Time to Respond to FDA <input type="checkbox"/> Request Meeting <input type="checkbox"/> Request Hearing			
<input type="checkbox"/> Report submission: <input type="checkbox"/> Current Investigator <input type="checkbox"/> Annual Progress Report <input type="checkbox"/> Site Waiver Report <input type="checkbox"/> Final					
<input type="checkbox"/> Other Reason ( <i>specify</i> )					

SECTION D3			REASON FOR SUBMISSION - 510(k)		
<input type="checkbox"/> New Device	<input type="checkbox"/> Additional or Expanded Indications	<input type="checkbox"/> Change in Technology			
<input checked="" type="checkbox"/> Other Reason ( <i>specify</i> ) change in sterilization method					

**SECTION E ADDITIONAL INFORMATION ON 510(K) SUBMISSIONS**

Product codes of devices to which substantial equivalence is claimed				Summary of, or statement concerning, safety and effectiveness information	
1	GEX	2		3	
5		6		7	
				<input checked="" type="checkbox"/> 510 (k) summary attached <input type="checkbox"/> 510 (k) statement	

Information on devices to which substantial equivalence is claimed (if known)

	510(k) Number		Trade or Proprietary or Model Name		Manufacturer
1	K100384	1	<b>FiberLase CO2 Laser WaveGuide</b>	1	<b>Lumenis</b>
2					
3					
4					
5		5		5	
6		6		6	

**SECTION F PRODUCT INFORMATION - APPLICATION TO ALL APPLICATIONS**

Common or usual name or classification  
**FiberLase Endure CO<sub>2</sub> Fiber**

	Trade or Proprietary or Model Name for This Device		Model Number
1	<b>FiberLase Endure CO<sub>2</sub> Fiber</b>	1	
2		2	
3		3	
4		4	

FDA document numbers of all prior related submissions (regardless of outcome)

1	2	3	4	5	6
7	8	9	10	11	12

Data Included in Submission  
 Laboratory Testing       Animal Trials       Human Trials

**SECTION G PRODUCT CLASSIFICATION - APPLICATION TO ALL APPLICATIONS**

Product Code GEX	C.F.R. Section (if applicable) 21 CFR 878.4810	Device Class <input type="checkbox"/> Class I <input checked="" type="checkbox"/> Class II <input type="checkbox"/> Class III <input type="checkbox"/> Unclassified
Classification Panel General and Plastic Surgery		

Indications (from labeling)  
 The Lumenis FiberLase ENDURE CO<sub>2</sub> fiber is intended for use in surgical procedures requiring ablation, vaporization, excision, incision and coagulation of soft tissue. The FiberLase CO<sub>2</sub> fiber is indicated for use in open surgical procedures such as ENT surgery and laparoscopy and endoscopic procedures

Note: Submission of this information does not affect the need to submit a 2891 or 2891a Device Establishment Registration form.	FDA Document Number <i>(if known)</i>
---	---------------------------------------

**SECTION H MANUFACTURING / PACKAGING / STERILIZATION SITES RELATING TO A SUBMISSION**

<input checked="" type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete	FDA Establishment Registration Number NA	<input checked="" type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Manufacturer	<input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager / Relabeler
Company / Institution Name Lumenis Ltd.		Establishment Registration Number NA	
Division Name <i>(if applicable)</i>		Phone Number <i>(including area code)</i> (972)4-959 9000	
Street Address 13 Hayetzira Street		FAX Number <i>(including area code)</i> (9724-959 9050)	
City Yokneam	State / Province	ZIP/Postal Code 20692	Country Israel
Contact Name Yoram levy	Contact Title Lumenis QA/RA Consultant	Contact E-mail Address Yoram@qsitemed.com	

<input checked="" type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete	FDA Establishment Registration Number (b) (4)	<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Manufacturer	<input checked="" type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager / Relabeler
--	--	---	--

(b) (4)

<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete	FDA Establishment Registration Number	<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Manufacturer	<input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager / Relabeler
Company / Institution Name		Establishment Registration Number	
Division Name <i>(if applicable)</i>		Phone Number <i>(including area code)</i> (    )	
Street Address		FAX Number <i>(including area code)</i> (    )	
City	State / Province	ZIP/Postal Code	Country
Contact Name	Contact Title	Contact E-mail Address	

**SECTION I****UTILIZATION OF STANDARDS**

Note: Complete this section if your application or submission cites standards or includes a "Declaration of Conformity to a Recognized Standard" statement.

	Standards No.	Standards Organization	Standards Title	Version	Date
1	ANSI/AAMI/ISO 11135-1	ANSI/AAMI/ISO	Sterilization of health care products — Ethylene oxide	NA	2007
2	AAMI TIR30:2001	AAMI	Product adoption and process equivalency for ethylene oxide sterilization	NA	2001
	Standards No.	Standards Organization	Standards Title	Version	Date
3	IEC 60601-2-22	IEC	Medical Electrical Equipment - Part 2-22: Particular Requirements for the Safety of Diagnostic and Therapeutic Laser Equipment.	NA	1995 amended on 2007
	Standards No.	Standards Organization	Standards Title	Version	Date
4	IEC 60825-1	IEC	Safety of Laser Products - Part 1: Equipment Classification, Requirements and User's Guide.	NA	2007
5	ANSI/AAMI/ISO 17665-1:2006	ANSI/AAMI/ISO	Sterilization of health care products — Moist Heat.	NA	2006
6	ISO 14971-1	ISO	Risk management for medical devices	NA	2007

Please include any additional standards to be cited on a separate page.

**Public reporting burden for this collection of information** is estimated to average 0.5 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to:

Food and Drug Administration  
CDRH (HFZ-342)  
9200 Corporate Blvd.  
Rockville, MD 20850

*An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control*

Department of Health and Human Services  
Food and Drug Administration

**STANDARDS DATA REPORT FOR 510(K)S**  
*(To be filled in by applicant)*

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional       Special       Abbreviated

STANDARD TITLE<sup>1</sup>

AANSI/AAMI/ISO 17665-1:2006 Sterilization of health care products - Moist heat - Part 1 Requirements for the development, validation and routine control of a sterilization process for medical devices

**Please answer the following questions**

Yes      No

Is this standard recognized by FDA<sup>2</sup>? .....      

FDA Recognition number<sup>3</sup> ..... # 0214-261

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report<sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of the standard? .....         
If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
If yes, were deviations in accordance with the FDA supplemental information sheet (SIS)<sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
If yes, report these exclusions in the summary report table.

Is there an FDA guidance<sup>6</sup> that is associated with this standard? .....         
If yes, was the guidance document followed in preparation of this 510k? .....      

Title of guidance: \_\_\_\_\_

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]  
<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)  
<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>  
<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>.

<sup>6</sup> The online search of CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

25

**EXTENT OF STANDARD CONFORMANCE  
SUMMARY REPORT TABLE**

STANDARD TITLE  
AANSI/AAMI/ISO 17665-1:2006 STERILIZATION OF HEALTH CARE PRODUCTS - MOIST HEAT - PART 1 REQUIREMENTS FOR THE DEVELOPMENT, VALIDATION AND ROUTINE CONTROL OF A STERILIZATION PROCESS FOR MEDICAL DEVICES

**CONFORMANCE WITH STANDARD SECTIONS\***

SECTION NUMBER	SECTION TITLE N/A There were no deviations or alternate test options used	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
----------------	--	--

TYPE OF DEVIATION OR OPTION SELECTED<sup>L</sup>

DESCRIPTION

JUSTIFICATION

SECTION NUMBER	SECTION TITLE	CONFORMANCE? <input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
----------------	---------------	---

TYPE OF DEVIATION OR OPTION SELECTED<sup>L</sup>

DESCRIPTION

JUSTIFICATION

SECTION NUMBER	SECTION TITLE	CONFORMANCE? <input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
----------------	---------------	---

TYPE OF DEVIATION OR OPTION SELECTED<sup>L</sup>

DESCRIPTION

JUSTIFICATION

\* For completeness list all sections of the standard and indicate whether conformance is met. If a section is not applicable (N/A) an explanation is needed under "justification." Some standards include options, so similar to deviations, the option chosen needs to be described and adequately justified as appropriate for the subject device. Explanation of all deviations or description of options selected when following a standard is required under "type of deviation or option selected," "description" and "justification" on the report. More than one page may be necessary.

<sup>L</sup> Types of deviations can include an exclusion of a section in the standard, a deviation brought out by the FDA supplemental information sheet (SIS), a deviation to adapt the standard to the device, or any adaptation of a section.

**Paperwork Reduction Act Statement**

Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to:

Center for Devices and Radiological Health  
1350 Piccard Drive  
Rockville, MD 20850

*An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number.*

26





## **Section C – 510k Cover Letter**



K-7

R130164

January 21, 2013  
Food and Drug Administration  
Center for Medical Devices and Radiological Health  
Office of Device Evaluation  
10903 New Hampshire Avenue  
White Oak, Bldg 66, Room 2534  
Silver Spring MD 20990  
USA

FDA CDRH DMC  
JAN 23 2013  
Received

**Re: Special Premarket Notification 510(k) Notification  
For *FiberLase Endure CO2 Fiber***

Dear Madam/Sir:

In accordance with section 510(k) of the Federal Food, Drug, and Cosmetic Act, and with the Code of Federal Regulations, 21 CFR 807, Subpart E, Lumenis Ltd. ("Lumenis" or "the Company"), hereby submits a special 510(k) Premarket Notification for *FiberLase Endure CO2 Fiber*.

For the convenience of reviewers, the Agency's "Screening Checklist for All Premarket Notifications" is included in this submission. The checklist identifies the sections within this submission where the required information is located.

As required by MDUFMA of 2002, a copy of the Medical Device User Fee Cover Sheet (Payment Identification Number (b) (4) ) is provided in Section A of this application.

Presented below is a summary of applicable information related to this submission. Detailed information is provided in Sections 1 to 16 of the submission.

**Applicant Name:** Lumenis Ltd.  
13 Hayetzira Street  
Yokneam Industrial Park

C-2

*FiberLase Endure CO2 Fiber - 510k Notification*

28

Detailed information is provided in Sections 1 to 16 of the submission.

**Applicant Name:** Lumenis Ltd.  
13 Hayetzira Street  
Yokneam Industrial Park  
Yokneam 20692 Israel  
Tel. (972)4-959-9000  
Fax: (972)4-959-9050

**Contact Person:** Yoram Levy, Qsite  
31 Haavoda St.  
Binyamina, ISRAEL  
Tel: (972)4-638-8837  
Fax: (972)4-638-0510

**Device Trade Name:** *FiberLase Endure CO2 Fiber*

**Device Type:** CO2 Laser fiber

**Classification:** **Regulation Name:** Laser surgical instrument for use in general and plastic surgery and in dermatology  
**Product Code:** GEX  
**Regulation No:** 21 CFR 878.4810  
**Class:** II  
**Classification Panel:** General & Plastic Surgery

**Indications for Use:** The *FiberLase Endure CO2 Fiber* is intended for use in surgical procedures requiring ablation, vaporization, excision, incision and coagulation of soft tissue. The *FiberLase Endure CO2 Fiber* is indicated for use in open surgical procedures such as ENT surgery and laparoscopy and endoscopic procedures.

**Predicate device:** **FiberLase CO2 Laser WaveGuide (K100384)**  
**Regulation Number:** 21 CFR 878.4810  
**Regulation Name:** Laser surgical instrument for use in  
general and plastic surgery and in dermatology  
**Regulatory Class:** II  
**Product Code:** GEX

**Type of submission:** Special 510k

**Design and Use of the Device:**

Question	YES	NO
Is the device intended for prescription use (21 CFR 801 Subpart D)? <sup>A</sup>	X	
Is the device intended for over-the-counter use (21 CFR 807 Subpart C)?		X
Does the device contain components derived from a tissue or other biologic source?		X
Is the device provided sterile?	X	
Is the device intended for single use?		X
Is the device a reprocessed single use device?		X
If yes, does this device type require reprocessed validation data?		N/A
Does the device contain a drug?		X
Does the device contain a biologic?		X
Does the device use software?		X
Does the submission include clinical information?		X
Is the device implanted?		X

**Confidentiality:** Lumenis considers its intent to market the *FiberLase Endure CO2 Fiber* to be confidential commercial information. The Company has not disclosed its intent to market this device to anyone except its employees, others with a financial interest in

the Company, its advertising and law firms, and its consultants. The Company, therefore, requests that FDA not disclose the existence of this application until such time as final action on the submission is taken. In addition, some of the material in this application may be trade secret or confidential commercial or financial information within the meaning of 21 CFR § 20.61 and therefore not disc losable under the Freedom of Information Act, even after the existence of the application becomes public. We ask that FDA consult with the Company as provided in 21 CFR § 20.45 before making any part of this submission publicly available.

We trust that the information contained in this 510(k) Notification will be sufficient to enable FDA to find that the *FiberLase Endure CO2 Fiber* is substantially equivalent to its predicate device.

Please direct any correspondence concerning this submission to Yoram Levy of QSite, Lumenis regulatory consultant, at [Yoram@qsitemed.com](mailto:Yoram@qsitemed.com), Tel (972)4-638-8837 or Fax (972)4-638-0510. Upon a finding of substantial equivalence, please send me a copy of the signed substantial equivalence letter by facsimile to (972)4-638-0510 or by Email.

Sincerely,

  
Yoram Levy  
Lumenis Ltd. QA/RA Consultant

CC: Assaf Gelstein, Lumenis Ltd. Platform Manager

\* \* \* COMMUNICATION RESULT REPORT ( MAR. 27. 2013 10:32AM ) \* \* \*

FAX HEADER 1:  
FAX HEADER 2:TRANSMITTED/STORED : MAR. 27. 2013 10:25AM  
E MODE OPTION

ADDRESS

RESULT

PAGE

4037 MEMORY TX

972 4 638 0510

OK

3/3

REASON FOR ERROR  
E-1) HANG UP OR LINE FAIL  
E-3) NO ANSWERE-2) BUSY  
E-4) NO FACSIMILE CONNECTION

DEPARTMENT OF HEALTH &amp; HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002Lumenis, Limited  
% Qsite  
Mr. Yoram Levy  
General Manager  
31 Haavoda Street  
Binyamina, Israel 30500.

March 24, 2012

Re: K130164

Trade/Device Name: FiberLase Endure CO2 Fiber  
Regulation Number: 21 CFR 878.4810  
Regulation Name: Laser surgical instrument for use in general and plastic surgery and in dermatology  
Regulatory Class: Class II  
Product Code: GEX  
Dated: March 06, 2013  
Received: March 11, 2013

Dear Mr. Levy:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set



Food and Drug Administration  
Office of Device Evaluation &  
Office of In Vitro Diagnostics

**COVER SHEET MEMORANDUM**

**\*\*\*NOTE: This form is OPTIONAL for holds, and REQUIRED for final decisions\*\*\***

**From:** Reviewer Name Xin Tan  
**Subject:** 510(k) Number K130164/S001  
**To:** The Record

Please list CTS decision code SE

- Refused to accept (Note: this is considered the first review cycle, See Screening Checklist [http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0\\_5631/Screening%20Checklist%207%202%2007.doc](http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0_5631/Screening%20Checklist%207%202%2007.doc))
- Hold (Additional Information or Telephone Hold).
- Final Decision (SE, SE with Limitations, NSE (select code below), Withdrawn, etc.).

Not Substantially Equivalent (NSE) Codes

- NO NSE for lack of predicate
- NI NSE for new intended use
- NQ NSE for new technology that raises new questions of safety and effectiveness
- NU NSE for new intended use AND new technology raising new questions of safety and effectiveness
- NP NSE for lack of performance data
- NS NSE no response
- NL NSE for lack of performance data AND no response
- NM NSE pre-amendment device call for PMAs (515i)
- NC NSE post-amendment device requires PMAs
- NH NSE for new molecular entity requires PMA
- TR NSE for transitional device

Please complete the following for a final clearance decision (i.e., SE, SE with Limitations, etc.):		YES	NO
Indications for Use Page	<i>Attach IFU</i>	X	
510(k) Summary /510(k) Statement	<i>Attach Summary</i>	X	
Truthful and Accurate Statement.	<i>Must be present for a Final Decision</i>	X	
Is the device Class III? If yes, does firm include Class III Summary?	<i>Must be present for a Final Decision</i>		X
Does firm reference standards? (If yes, please attach form from <a href="http://www.fda.gov/opacom/morechoices/fdaforms/FDA-3654.pdf">http://www.fda.gov/opacom/morechoices/fdaforms/FDA-3654.pdf</a> )		X	
Is this a combination product? (Please specify category _____, see <a href="http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0_413b/COMBINATION%20PRODUCT%20ALGORITHM%20(REVISED%203-12-03).DOC">http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0_413b/COMBINATION%20PRODUCT%20ALGORITHM%20(REVISED%203-12-03).DOC</a> )			X
Is this a reprocessed single use device? (Guidance for Industry and FDA Staff – MDUFMA - Validation Data in 510(k)s for Reprocessed Single-Use Medical Devices, <a href="http://www.fda.gov/cdrh/ode/guidance/1216.html">http://www.fda.gov/cdrh/ode/guidance/1216.html</a> )			X
Is this device intended for pediatric use only?			X
Is this a prescription device? (If both prescription & OTC, check both boxes.)		X	
Did the application include a completed FORM FDA 3674, Certification with Requirements of ClinicalTrials.gov Data Bank?			X
Is clinical data necessary to support the review of this 510(k)?			X

Rev. 9/20/12 – added digital concurrence table



For United States-based clinical studies only: Did the application include a completed FORM FDA 3674, <i>Certification with Requirements of ClinicalTrials.gov Data Bank?</i> (If study was conducted in the United States, and FORM FDA 3674 was not included or incomplete, then applicant must be contacted to obtain completed form.)			
Does this device include an Animal Tissue Source?			X
All Pediatric Patients age <=21			X
Neonate/Newborn (Birth to 28 days)			X
Infant (29 days -< 2 years old)			X
Child (2 years -< 12 years old)			X
Adolescent (12 years -< 18 years old)			X
Transitional Adolescent A (18 - <21 years old) Special considerations are being given to this group, different from adults age ≥ 21 (different device design or testing, different protocol procedures, etc.)			X
Transitional Adolescent B (18 -<= 21; No special considerations compared to adults => 21 years old)			X
Nanotechnology			X
Mobile Application			X
MR Conditional			X
Device Contains Battery			X
Companion Diagnostic			X
Is this device subject to the Tracking Regulation? (Medical Device Tracking Guidance, <a href="http://www.fda.gov/cdrh/comp/guidance/169.html">http://www.fda.gov/cdrh/comp/guidance/169.html</a> )	Contact OC.		X

**Regulation Number**                      **Class\***                      **Product Code**  
 21 CFR 878.4810                      Class II                      GEX

(\*If unclassified, see 510(k) Staff)

**Additional Product Codes:** \_\_\_\_\_

Digital Signature Concurrence Table	
Reviewer Sign-Off	Xin Tan -S 2013.03.22 12:20:10 -04'00'
Branch Chief Sign-Off	Neil R Ogden 2013.03.22 16:14:20 -04'00'
Division Sign-Off	Peter D. Rumm -S 2013.03.25 22:06:32 -04'00'

**Tan, Xin**

---

**From:** Yoram Levy <yoram@qsitemed.com>  
**Sent:** Monday, March 11, 2013 3:19 AM  
**To:** Tan, Xin  
**Subject:** RE: Your 510(k) Submission K130164: Additional Information Required

Dear Sofia,

I organized your questions in the two last emails in a response format so I will be able to address all of them:

(b)(4)



(b)(4)



Please email any farther questions that you have.

Best Regards,

*Yoram*

Yoram Levy, Qsite



*Global Quality & Regulatory*

Tel (972)4-638-8837, Fax (972)4-638-0510  
Cell (972)52-279-2871



Food and Drug Administration  
Office of Device Evaluation &  
Office of In Vitro Diagnostics

**COVER SHEET MEMORANDUM**

**From:** Reviewer Name Xin TAN  
**Subject:** 510(k) Number K130164  
**To:** The Record

Please list CTS decision code AI

- Refused to accept (Note: this is considered the first review cycle, See Screening Checklist [http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0\\_5631/Screening%20Checklist%207%20202%2007.doc](http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0_5631/Screening%20Checklist%207%20202%2007.doc))
- Hold (Additional Information or Telephone Hold).
- Final Decision (SE, SE with Limitations, NSE (select code below), Withdrawn, etc.).

Not Substantially Equivalent (NSE) Codes

- NO NSE for lack of predicate
- NI NSE for new intended use
- NQ NSE for new technology that raises new questions of safety and effectiveness
- NU NSE for new intended use AND new technology raising new questions of safety and effectiveness
- NP NSE for lack of performance data
- NS NSE no response
- NL NSE for lack of performance data AND no response
- NM NSE pre-amendment device call for PMAs (515i)
- NC NSE post-amendment device requires PMAs
- NH NSE for new molecular entity requires PMA
- TR NSE for transitional device

Please complete the following for a final clearance decision (i.e., SE, SE with Limitations, etc.):		YES	NO
Indications for Use Page	Attach IFU		
510(k) Summary /510(k) Statement	Attach Summary		
Truthful and Accurate Statement.	Must be present for a Final Decision		
Is the device Class III?			
If yes, does firm include Class III Summary?	Must be present for a Final Decision		
Does firm reference standards? (If yes, please attach form from <a href="http://www.fda.gov/opacom/morechoices/fdaforms/FDA-3654.pdf">http://www.fda.gov/opacom/morechoices/fdaforms/FDA-3654.pdf</a> )			
Is this a combination product? (Please specify category _____, see <a href="http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0_413b/CO-MBINATION%20PRODUCT%20ALGORITHM%20(REVISED%203-12-03).DOC">http://eroom.fda.gov/eRoomReq/Files/CDRH3/CDRHPremarketNotification510kProgram/0_413b/CO-MBINATION%20PRODUCT%20ALGORITHM%20(REVISED%203-12-03).DOC</a> )			
Is this a reprocessed single use device? (Guidance for Industry and FDA Staff – MDUFMA - Validation Data in 510(k)s for Reprocessed Single-Use Medical Devices, <a href="http://www.fda.gov/cdrh/ode/guidance/1216.html">http://www.fda.gov/cdrh/ode/guidance/1216.html</a> )			
Is this device intended for pediatric use only?			
Is this a prescription device? (If both prescription & OTC, check both boxes.)			
Did the application include a completed FORM FDA 3674, Certification with Requirements of ClinicalTrials.gov Data Bank?			
Is clinical data necessary to support the review of this 510(k)?			
For United States-based clinical studies only: Did the application include a completed FORM FDA 3674, Certification with Requirements of ClinicalTrials.gov Data Bank? (If study was			

conducted in the United States, and FORM FDA 3674 was not included or incomplete, then applicant must be contacted to obtain completed form.)

Does this device include an Animal Tissue Source?

All Pediatric Patients age <= 21

Neonate/Newborn (Birth to 28 days)

Infant (29 days - < 2 years old)

Child (2 years - < 12 years old)

Adolescent (12 years - < 18 years old)

Transitional Adolescent A (18 - < 21 years old) Special considerations are being given to this group, different from adults age >= 21 (different device design or testing, different protocol procedures, etc.)

Transitional Adolescent B (18 - <= 21; No special considerations compared to adults => 21 years old)

Nanotechnology

Is this device subject to the Tracking Regulation? (Medical Device Tracking Guidance, <http://www.fda.gov/cdrh/comp/guidance/169.html>) Contact OC.

Regulation Number

21 CFR 878.4210

Class\*

Class II

Product Code

GEX

(\*If unclassified, see 510(k) Staff)

Additional Product Codes:

Review:

*[Signature]*

(Branch Chief)

G52B1

(Branch Code)

2/27/13


(Date)

Final Review:

(Division Director)

(Date)

### 510(k) "SUBSTANTIAL EQUIVALENCE" DECISION-MAKING PROCESS



\* 510(k) Submissions compare new devices to marketed devices. FDA requests additional information if the relationship between marketed and "predicate" (pre-Amendments or reclassified post-Amendments) devices is unclear.

\*\* This decision is normally based on descriptive information alone, but limited testing information is sometimes required.

\*\*\* Data maybe in the 510(k), other 510(k)s, the Center's classification files, or the literature.



**DEPARTMENT OF HEALTH AND HUMAN SERVICES      M E M O R A N D U M**

Food and Drug Administration  
Office of Device Evaluation  
WO Building 66  
Silver Spring, MD 20993

**Premarket Notification [510(k)] Review  
Special  
K130164**

**Date:** February 15, 2013  
**To:** The Record  
**From:** (Sofia) Xin Tan, Ph.D., PHYSICIST  
**Office:** ODE  
**Division:** DSD/GSDB1

**510(k) Holder:** Lumenis Ltd. (Israel)  
**Device Name:** *FiberLase Endure CO<sub>2</sub> Fiber*

**Contact:** Yoram Levy  
Qsite  
31 Haavoda Street  
Binyamina, ISRAEL  
Phone: (972) 4-638-8837  
Fax: (972) 4-638-0510  
Email: yoram@qsitemed.com

**I. Purpose & Submission Summary**

The 510(k) holder (Lumenis Ltd.) submitted the Special Premarket Notification (K130164) to request the U.S. marketing clearance for the proposed device (*FiberLase Endure CO<sub>2</sub> Fiber*). The company's own FDA-cleared device (K100384 *Lumenis FiberLase CO<sub>2</sub> Laser WaveGuide*) is identified as the predicate. The sponsor claimed Substantial Equivalence based on the same intended use and clinical indications, same fundamental scientific technology and materials, with the only modification in changing the fiber from single to multi use by addition of an autoclave sterilization method for reprocessing the fiber.

On 14 February 2013, the sponsor sent via email a revised 510(k) Cover Letter with the added statement that there were no prior submissions for the proposed device.

**II. Device Description & Indications for Use**

	Yes	No	N/A
Is the device life-supporting or life sustaining?		X	

Is the device an implant (implanted longer than 30 days)?		X	
Is the device intended for prescription use (21 CFR 801 Subpart D)?	X		
Is the device intended for over-the-counter use (21 CFR 807 Subpart C)?		X	
Does the device use software?		X	
Is the device sterile?	X		
Is the device reusable (not reprocessed single use)?	X		
Are "cleaning" instructions included for the end user? Cleaning and autoclave sterilization methods were provided for reprocessing.	X		

(b) (4)



Submission Number	Device Name	Sponsor	Indications for Use
<b>Subject</b> K130164	<i>FiberLase Endure CO<sub>2</sub> Fiber</i>	Lumenis Ltd.	(b) (4)
<b>Predicate</b> K100384	<i>FiberLase CO<sub>2</sub> Laser WaveGuide</i>	Lumenis Ltd.	



**III. Predicate Device Comparison**

(b) (4)



**IV. Labeling**

The sponsor provided comprehensive User Manuals for the *FiberLase Endure CO<sub>2</sub> Fiber*, the Handpieces, and the Renewal Kit. The FiberLase Handpieces are designed to facilitate the delivery of the fiber to the tissue treatment site. The FiberLase Fiber Renewal Kit is designed to repair the fiber tip if it becomes damaged during a surgical CO<sub>2</sub> laser procedure. The User Manuals also contain instructions for cleaning and sterilization. The labeling section is adequate.

**V. Sterilization/Shelf Life**

(b) (4)



**VI. Biocompatibility**

(b) (4)



**VII. Software**


Not applicable. The subject device does not use any software.

**VIII. Electromagnetic Compatibility, Electrical and Mechanical Safety**

Not applicable. The device does not contain Electrical and Electromagnetic Systems.

**IX. Performance Testing - Bench**

(b) (4)



**X. Performance Testing - Animal**

Not applicable. This submission does not include animal testing or data.

**XI. Performance Testing - Clinical**

Not applicable. This submission does not include human clinical testing or data.

**XII. Substantial Equivalence Discussion**

	Yes	No	
1. Is Product A Device	X		If NO = Stop, see 510(k) staff
2. Is Device Subject To 510(k)?	X		If NO = Stop, see 510(k) staff
3. Same Indication Statement?	X		If YES = Go To 5
4. Do Differences Alter The Effect Or Raise New Issues of Safety Or Effectiveness?			If YES = Stop NSE
5. Same Technological Characteristics?	X		If YES = Go To 7
6. Could The New Characteristics Affect Safety Or Effectiveness?			If YES = Go To 8
7. Descriptive Characteristics Precise Enough? Need performance data to evaluate substantial equivalence after the addition of an autoclave (steam) sterilization procedure to reprocess the device for multi-use.		X	If NO = Go To 10 If YES = Stop SE
8. New Types Of Safety Or Effectiveness Questions?			If YES = Stop NSE
9. Accepted Scientific Methods Exist?			If NO = Stop NSE
10. Performance Data Available?	X		If NO = Request Data

18

II. Data Demonstrate Equivalence? Questions pertaining to performance testing need to be addressed before determination on Substantial Equivalence.	?	Final Decision:
--	---	-----------------

Note: Please complete the above table and answer the corresponding questions. "Yes" responses to questions 4, 6, 8, and 11, and every "no" response requires an explanation.

**XIII. Deficiencies**

(b) (4)



**XIV. Recommendation**

(b) (4)

**Regulation Number:** 21 CFR 878.4810

**Regulation Name:** Laser surgical instrument for use in general and plastic surgery and in dermatology

**Device Class:** Class II

**Product Code:** GEX

Xin Tan - S  2013.02.26  
10:27:36 -05'00'

---

Xin Tan  
**Lead Reviewer**

Neil R Ogden  2013.02.26  
15:24:24 -05'00'

---

Neil Ogden  
**Branch Chief, ODE/DSD/GSDB1**

re

**Ogden, Neil**

---

**From:** Tan, Xin  
**Sent:** Tuesday, February 26, 2013 3:36 PM  
**To:** 'Yoram Levy'  
**Cc:** Ogden, Neil  
**Subject:** Your 510(k) Submission K130164: Additional Information Required

Yoram,

(b) (4)



Performance Testing

(b) (4)



Sterilization

(b) (4)



Please submit your full response to the above deficiencies via email at your earliest convenience. Let me know if you have any questions regarding the content of this email.

2/

Sincerely,

Sofia

---

**XIN TAN, Ph.D. PHYSICIST**

U.S. Food and Drug Administration

CDRH/OSEL/DP

10903 New Hampshire Avenue

White Oak Building 62, Room 1110

Silver Spring, Maryland 20993-0002

Phone: 301-796-2719

Fax: 301-796-9927

Email: [xin.tan@fda.hhs.gov](mailto:xin.tan@fda.hhs.gov)

K130164/S001



FDA CDRH DMC

MAR 11 2013

Received

March 6, 2013  
Food and Drug Administration  
Center for Medical Devices and Radiological Health  
Office of Device Evaluation  
10903 New Hampshire Avenue  
White Oak, Bldg 66, Room 2534  
Silver Spring MD 20990  
USA

**Re: Response for the Lumenis *FiberLase Endure* eCopy Cover Letter**

Dear Madam/Sir:

In accordance with the eCopy Program for Medical Device Submissions Lumenis, ("the Company") is submitting the enclosed response for 510(k) Notification (21 CFR 807.90(e)) eCopy that escorts the printed hard copy of the same response.

**The eCopy is an exact duplicate of the paper copy.**

Presented below is a summary of applicable information related to this submission. Detailed information is provided in Sections 1 to 16 of the submission.

**Applicant Name:**

Lumenis Ltd.  
13 Hayetzira Street  
Yokneam Industrial Park  
Yokneam 20692 Israel  
Tel. (972)4-959 9000  
Fax: (972)4-959 9050

**Contact Person:**

Yoram Levy  
Haavoda 31,  
Binyamina 30500  
Israel  
Phone: (972)4-638-8837  
Fax: (972)4-638-0510  
E-mail: [Yoram@qsitemed.com](mailto:Yoram@qsitemed.com)

**Trade Name:**

**Lumenis *FiberLase Endure***


Please direct any correspondence concerning this submission to Yoram Levy of QSite, Lumenis regulatory consultant, at [Yoram@qsitemed.com](mailto:Yoram@qsitemed.com), Tel (972)4-638-8837 or Fax (972)4-638-0510. Upon a finding of substantial equivalence, please send me a copy of the signed substantial equivalence letter by facsimile to (972)4-638-0510 or by Email.

***Lumenis ResurFX 510k Notification Cover Letter***

52



Sincerely,

  
Yoram Levy  
Lumenis Ltd. QA/RA Consultant

***Lumenis ResurFX 510k Notification Cover Letter***





Date: March 3, 2013

To: Sofia XIN TAN, Ph.D. PHYSICIST  
U.S. Food and Drug Administration  
CDRH/OSEL/DP  
10903 New Hampshire Avenue  
White Oak Building 62, Room 1110  
Silver Spring, Maryland 20993-0002  
USA

**Subject: Response to K130164 Fiberlase Endure CO2 Fiber Request for Additional Information**

Dear Dr. XIN TAN

(b) (4)





(b) (4)

A large, solid grey rectangular redaction box covers the majority of the page's content, starting below the logo and extending nearly to the bottom of the page.

(b) (4)

A horizontal, solid grey rectangular redaction bar is located near the bottom of the page, below the large redaction box.



(b) (4)

A large, solid grey rectangular redaction box covers the majority of the page's content, starting below the LUMENIS logo and extending nearly to the bottom of the page. The text '(b) (4)' is positioned at the top left corner of this redacted area.



(b) (4)

A large, solid grey rectangular box covers the majority of the page, indicating that the content has been redacted under FOIA exemption (b)(4).



(b) (4)

A large, solid grey rectangular redaction box covers the majority of the page, obscuring all text and graphics that would otherwise be present.



(b) (4)



Sincerely,

A handwritten signature in black ink, appearing to read 'Yoram Levy'.

Yoram levy

[yoram@qsitemed.com](mailto:yoram@qsitemed.com)

(972)4-638-8837 Cell (972)52-279-2871

Fax (972)4-638-0510



## Attachment Table of Contents

Attachment No. 1: (b) (4) Waveguide Data Sheet ..... 8

Attachment No. 2: ANSI/AAMI/ISO 17665-1:2006 Sterilization of Health Care Products – Moist Heat Standards Data Report Form (FDA 3654) ..... 9



**Attachment No. 1:** (b) (4) **Waveguide Data Sheet**





**Attachment No. 2: ANSI/AAMI/ISO 17665-1:2006  
Sterilization of Health Care Products – Moist Heat  
Standards Data Report Form (FDA 3654)**



(b) (4) [Redacted]

Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(K)**  
*(To be filled in by applicant)*

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional                       Special                       Abbreviated

STANDARD TITLE <sup>1</sup>  
 AANSI/AAMI/ISO 17665-1:2006 Sterilization of health care products - Moist heat - Part 1 Requirements for the development, validation and routine control of a sterilization process for medical devices

**Please answer the following questions** Yes      No

Is this standard recognized by FDA <sup>2</sup>? .....      

FDA Recognition number <sup>3</sup> ..... # 0214-261

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report <sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
 If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
 If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of the standard? .....         
 If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
 If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) <sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
 If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
 If yes, report these exclusions in the summary report table.

Is there an FDA guidance <sup>6</sup> that is associated with this standard? .....         
 If yes, was the guidance document followed in preparation of this 510k? .....      

Title of guidance: \_\_\_\_\_

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]  
<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)  
<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>  
<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name address of the test laboratory or certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.  
<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>  
<sup>6</sup> The online search of CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

**EXTENT OF STANDARD CONFORMANCE  
SUMMARY REPORT TABLE**

STANDARD TITLE  
AANSI/AAMI/ISO 17665-1:2006 STERILIZATION OF HEALTH CARE PRODUCTS - MOIST HEAT - PART 1 REQUIREMENTS FOR THE DEVELOPMENT, VALIDATION AND ROUTINE CONTROL OF A STERILIZATION PROCESS FOR MEDICAL DEVICES

**CONFORMANCE WITH STANDARD SECTIONS\***

SECTION NUMBER	SECTION TITLE N/A There were no deviations or alternate test options used	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
----------------	--	--

TYPE OF DEVIATION OR OPTION SELECTED<sup>□</sup>

DESCRIPTION

JUSTIFICATION

SECTION NUMBER	SECTION TITLE	CONFORMANCE? <input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
----------------	---------------	---

TYPE OF DEVIATION OR OPTION SELECTED<sup>□</sup>

DESCRIPTION

JUSTIFICATION

SECTION NUMBER	SECTION TITLE	CONFORMANCE? <input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
----------------	---------------	---

TYPE OF DEVIATION OR OPTION SELECTED<sup>□</sup>

DESCRIPTION

JUSTIFICATION

\* For completeness list all sections of the standard and indicate whether conformance is met. If a section is not applicable (N/A) an explanation is needed under "justification." Some standards include options, so similar to deviations, the option chosen needs to be described and adequately justified as appropriate for the subject device. Explanation of all deviations or description of options selected when following a standard is required under "type of deviation or option selected," "description" and "justification" on the report. More than one page may be necessary.

□ Types of deviations can include an exclusion of a section in the standard, a deviation brought out by the FDA supplemental information sheet (SIS), a deviation to adapt the standard to the device, or any adaptation of a section.

**Paperwork Reduction Act Statement**

Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to:

Center for Devices and Radiological Health  
1350 Piccard Drive  
Rockville, MD 20850

*An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number.*

**From:** Tan, Xin  
**Sent:** Friday, April 05, 2013 5:57 PM  
**To:** k130164@docs.fda.gov  
**Subject:** Email Exchange\_Additional Information and Response to Deficiencies

**From:** Yoram Levy [<mailto:yoram@qsitemed.com>]  
**Sent:** Sunday, March 03, 2013 12:37 PM  
**To:** Tan, Xin  
**Cc:** Ogden, Neil  
**Subject:** RE: Your 510(k) Submission K130164: Additional Information Required

Dear Dr. Xin Tan,

We responded to your request for additional information in the attached files. The attachments (include the form that was submitted again) are part of the "K130164 FiberLase Endure CO2 Fiber 510(k) Response".

Please email or call for any additional information that you may need for this 510k.

Best Regards,

*Yoram*  
Yoram Levy, Qsite  
**qsite headerN**



Tel (972)4-638-8837, Fax (972)4-638-0510  
Cell (972)52-279-2871

---

**From:** Tan, Xin [<mailto:Xin.Tan@fda.hhs.gov>]  
**Sent:** Tuesday, February 26, 2013 10:36 PM  
**To:** 'Yoram Levy'  
**Cc:** Ogden, Neil  
**Subject:** Your 510(k) Submission K130164: Additional Information Required

Yoram,

We have reviewed your 510(k) submission K130164 to market the device – FiberLase Endure CO2 Fiber. We cannot determine if the device is substantially equivalent to your legally marketed predicate device (FiberLase CO2 Laser WaveGuide) based solely on the information you provided. To complete the review of your submission, we require the following information:

Performance Testing

(b) (4) [Redacted]

(b) (4)

A large area of the document is redacted with grey bars. The redaction covers approximately 10 lines of text, starting from the first line below the header and ending just above the main body of text.

(b) (4)

A very large rectangular area of the document is redacted with a solid grey block, covering the majority of the page's content.

Please submit your full response to the above deficiencies via email at your earliest convenience. Let me know if you have any questions regarding the content of this email.

Sincerely,

Sofia

---

**XIN TAN, Ph.D.** PHYSICIST  
U.S. Food and Drug Administration  
CDRH/OSEL/DP  
10903 New Hampshire Avenue  
White Oak Building 62, Room 1110  
Silver Spring, Maryland 20993-0002  
Phone: 301-796-2719  
Fax: 301-796-9927  
Email: [xin.tan@fda.hhs.gov](mailto:xin.tan@fda.hhs.gov)

**From:** Yoram Levy [mailto:[yoram@qsitemed.com](mailto:yoram@qsitemed.com)]  
**Sent:** Monday, March 11, 2013 8:55 PM  
**To:** Tan, Xin  
**Subject:** RE: Your 510(k) Submission K130164: Additional Information Required

Yes we did.

(b) (4)

Best Regards,

*Yoram*

Yoram Levy, Qsite

qsite headerN



Tel (972)4-638-8837, Fax (972)4-638-0510  
Cell (972)52-279-2871

---

**From:** Tan, Xin [mailto:[Xin.Tan@fda.hhs.gov](mailto:Xin.Tan@fda.hhs.gov)]  
**Sent:** Monday, March 11, 2013 11:47 PM  
**To:** 'Yoram Levy'  
**Subject:** RE: Your 510(k) Submission K130164: Additional Information Required

Yoram,

(b) (4)

Sofia

---

**From:** Yoram Levy [mailto:[yoram@qsitemed.com](mailto:yoram@qsitemed.com)]  
**Sent:** Monday, March 11, 2013 3:19 AM  
**To:** Tan, Xin  
**Subject:** RE: Your 510(k) Submission K130164: Additional Information Required

Dear Sofia,

(b) (4)

(b) (4)



(b) (4)





Please email any farther questions that you have.

Best Regards,

*Yoram*  
Yoram Levy, Qsite  
**qsite headerN**



Tel (972)4-638-8837, Fax (972)4-638-0510  
Cell (972)52-279-2871

---

**From:** Tan, Xin [<mailto:Xin.Tan@fda.hhs.gov>]  
**Sent:** Saturday, March 09, 2013 12:54 AM  
**To:** 'Yoram Levy'  
**Subject:** RE: Your 510(k) Submission K130164: Additional Information Required

Yoram,

(b) (4)

Sofia

---

**From:** Yoram Levy [<mailto:yoram@qsitemed.com>]  
**Sent:** Friday, March 08, 2013 10:49 AM  
**To:** Tan, Xin  
**Subject:** RE: Your 510(k) Submission K130164: Additional Information Required

Dear Sofia,

(b) (4)

Best Regards,

*Yoram*  
Yoram Levy, Qsite  
**qsite headerN**



Tel (972)4-638-8837, Fax (972)4-638-0510  
Cell (972)52-279-2871

**From:** Tan, Xin [<mailto:Xin.Tan@fda.hhs.gov>]  
**Sent:** Friday, March 08, 2013 1:08 AM  
**To:** 'Yoram Levy'  
**Subject:** RE: Your 510(k) Submission K130164: Additional Information Required

Yoram,

(b) (4)  
[Redacted]

[Redacted]

Sofia

---

**From:** Yoram Levy [<mailto:yoram@qsite.com>]  
**Sent:** Sunday, March 03, 2013 12:37 PM  
**To:** Tan, Xin  
**Cc:** Ogden, Neil  
**Subject:** RE: Your 510(k) Submission K130164: Additional Information Required

Dear Dr. Xin Tan,

(b) (4)  
[Redacted]

Best Regards,

*Yoram*  
Yoram Levy, Qsite  


Tel (972)4-638-8837, Fax (972)4-638-0510  
Cell (972)52-279-2871

---

**From:** Tan, Xin [<mailto:Xin.Tan@fda.hhs.gov>]  
**Sent:** Tuesday, February 26, 2013 10:36 PM  
**To:** 'Yoram Levy'  
**Cc:** Ogden, Neil  
**Subject:** Your 510(k) Submission K130164: Additional Information Required

Yoram,

(b) (4)  
[Redacted]

(b) (4) [Redacted]

Performance Testing

(b) (4) [Redacted]

(b) (4) [Redacted]

Please submit your full response to the above deficiencies via email at your earliest convenience. Let me know if you have any questions regarding the content of this email.

Sincerely,

Sofia

---  
**XIN TAN, Ph.D.** PHYSICIST  
U.S. Food and Drug Administration  
CDRH/OSEL/DP  
10903 New Hampshire Avenue  
White Oak Building 62, Room 1110  
Silver Spring, Maryland 20993-0002  
Phone: 301-796-2719  
Fax: 301-796-9927  
Email: [xin.tan@fda.hhs.gov](mailto:xin.tan@fda.hhs.gov)

