510(k) Summary

1. Date of Summary

October 10, 2013

2. 510(k) Applicant

Broncus Medical, Inc. 1400 N. Shoreline Blvd, Suite A8 Mountain View, California 94043 Phone: (650) 428-1600 FAX: (650) 428-1542

OCT 1 5 2013

Contact Person:	Gary Kaplan
Phone:	(650) 428-1600
Fax:	(650) 428-1542
e-mail:	gkaplan@broncus.com

3. Device Overview

Trade Name:	LungPoint [™] Tools (LungPoint Sheath and LungPoint Dilation Balloon)
Common Name:	Sheath and Dilation Balloon
Classification Name:	Bronchoscope and Accessories 21 CFR 874.4680 Product Code EOQ

4. Predicate Device

The predicate devices identified are as follows:

Trade Name	510(k) Submitter	510(k) Number
CRE Pulmonary Balloon Dilatation Catheter	Boston Scientific Corporation	K023337, cleared to market on November 18, 2002
Olympus Guide Sheath	Olympus Medical Systems Corporation	K060243, cleared to market on June 23, 2006

5. Device Description

The LungPoint Tools are endoscopic tools used during bronchoscopy procedures. The LungPoint Sheath is designed to be used with a bronchoscope to provide a working

Broncus Medical, Inc.

Page 1 of 4

channel through which endoscopic tools, such as needles, dilation balloons, or other endoscopic devices may be introduced to the targeted lung tissue within the respiratory organs. The LungPoint Dilation Balloon is used to dilate tissue of the bronchial tree and may be inserted through the sheath or directly through the working channel of the bronchoscope. The bronchoscope is advanced to a predefined target following guidance of the LungPoint Software (previously cleared under 510(k) #K112051 to guide endoscopic tools or catheters in the lungs or to enable marker placement in soft lung tissue).

The materials used in the LungPoint Tools are commonly used medical grade materials and include platinum iridium markers.

Catalog Number	Working Length	Maximum Catheter Outer Diameter (OD)	Catheter Internal Diameter (ID)	Minimum accessory length for use through sheath
10007-1	900mm	2.65mm	2.0mm	965mm

The LungPoint Sheath has the following specifications:

The LungPoint Dilation Balloon has the following specifications:

Catalog Number	Balloon Size (OD x length)	Rated Balloon Pressure	Maximum Catheter Outer Diameter (OD)	Catheter Length/working length
10008-1	4mm x 6mm	20atm	<u> </u>	1430mm/975mm

6. Intended Use

The LungPoint Tools are endoscopic tools intended to be used with LungPoint Software guided bronchoscopes. The LungPoint Sheath is intended to be used as a working channel through which endoscopic tools may be introduced to targeted tissue. The LungPoint Dilation Balloon is intended to dilate tissue of the bronchial tree and may be inserted through the LungPoint Sheath or directly through the working channel of the bronchoscope. Not for pediatric use.

7. Comparison to Predicate Device

The LungPoint Tools are commonly used endoscopic tools with the same technological characteristics as the predicate devices. The indications for use of the LungPoint Tools for use with a bronchoscope guided by the LungPoint Software are all within the intended use of the predicate devices, which is to aid in reaching a target in the respiratory organ

either directly as is the case with the sheath or through dilating target tissue with the dilation balloon. The technological characteristics of the subject devices are the same as those of the predicate devices with the following exceptions

• LungPoint Balloon: balloon size and length, balloon burst pressure and balloon catheter length as outlined.

	CRE Pulmonary Balloon Dilatation Catheter (K023337)	LungPoint Dilation Balloon
Balloon Size and Inflation (balloon pressure)	8mm @ 3 ATM 9mm @ 5.5 ATM 10mm @ 9 ATM	4mm @ 10 ATM
Catheter Length (cm)	155	143
Balloon Length (cm)	3.0	0.6
Rated Burst Pressure (atmospheres)	9	20

• LungPoint Sheath: stylet, catheter ID/OD and catheter length as outlined.

	Olympus Guide Sheath (K060243)	LungPoint Sheath
Stylet Provided	No	Yes - used to enhance pushability and to prevent airway mucosa entering sheath
Catheter ID/OD	2.1/2.7 mm	2.0/2.6 mm
Catheter Length	900 mm	975 mm

None of these differences raise new questions of safety and effectiveness. Performance of the subject devices has been verified by use of accepted methods.

8. Performance Data

The design and safety of the LungPoint Tools were verified by performing functional and performance testing. All tests were designed to subject the sheath and dilation balloon to stresses that exceed those which would be encountered during clinical use. Testing included the following:

- Dimensional testing
- Joint/tensile test
- Simulated use
- Balloon fatigue/burst pressure
- Balloon deflation time
- Radiopacity.

All testing results met the pre-determined acceptance criteria that were established in the test protocols. Based on the testing the LungPoint Tools are as safe, as effective, and perform at least as safely and effectively as the predicate devices.

Broncus Medical, Inc.

Page 3 of 4

9. Safety and Effectiveness

The LungPoint Tools labeling contains instructions for use and any necessary cautions and warnings, to provide for safe and effective use of the devices. The biocompatibility assessment of all patient contacting materials was performed in accordance with ISO 10993, *Biological Evaluation of Medical Devices*. Specifically, cytotoxicity, sensitization, intracutaneous reactivity and systemic toxicity (acute) were tested. In addition, the devices are sterilized using e-beam sterilization.

10. Conclusion

Based on the testing the LungPoint Tools are as safe, as effective, and perform at least as safely and effectively as the predicate devices.

Broncus Medical, Inc.

Page 4 of 4

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016




DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration 10903 New Hampshire Avenue Document Control Center – WO66-G609 Silver Spring, MD 20993-002

October 15, 2013



Re: K131234

Trade/Device Name: LungPoint[™] Tools (LungPoint Sheath and LungPoint Dilation Balloon)
Regulation Number: 21 CFR 874.4680
Regulation Name: Bronchoscope (Flexible or Rigid) and Accessories
Regulatory Class: Class II
Product Code: EOQ
Dated: August 26, 2013
Received: August 27, 2013

Dear (b)(6)

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016

Page 2 - Mr. (b)(6)

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <u>http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm</u> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm.

Sincerely yours,

Eric A. Mann -S

for Malvina B. Eydelman, M.D. Director Division of Ophthalmic and Ear, Nose and Throat Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

Indications for Use	
510(k) Number (if known):	K131234
Device Name:	LungPoint Tools - LungPoint Dilation Balloon and LungPoint Sheath
Indications for Use:	The LungPoint Tools are endoscopic tools intended to be used with LungPoint Software guided bronchoscopes. The LungPoint Sheath is intended to be used as a working channel through which endoscopic tools may be introduced to targeted tissue. The LungPoint Dilation Balloon is intended to dilate tissue of the bronchial tree and may be inserted through the LungPoint Sheath or directly through the working channel of the bronchoscope. Not for pediatric use.

Prescription Use X (Part 21 CFR 801 Subpart D)

AND/OR


Over-The-Counter Use _____ (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)


Concurrence of CDRH, Office of Device Evaluation (ODE)

Page ____ of ____

Srinivas Nandkumar -S



Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhsigov or 301-796-8118



Records processed under FOIA Request # 2015-7253; Rele

Date: April 30, 2013

K131234 FDA CORH DMC MAY 01 2013

Receivre

U.S. Food and Drug Administration Center for Devices and Radiological Heath Document Mail Center – WO66-0609 10903 New Hampshire Avenue Silver Spring, MD 20993-0002

RE: Premarket Notification 510(k) Electronic Copy:

Broncus Medical, Inc. LungPoint Tools

To Whom It May Concern:

The electronic files (eCopy) are an exact duplicate of the paper copy included with this submission.

If you should have any questions regarding this submission please contact me at (b)(4) Third Party Information br email (b)(4) Third Party Information Please fax any correspondence regarding this submission to (b)(4) Third Party Information b(4) Third Party Information

Sincerely,

D)(6)

(b)(4) Third Party Information

(b)(4) Revision DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

(b)(4) Third Party Information

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016



Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

Date: April 29, 2013

U.S. Food and Drug Administration Center for Devices and Radiological Heath Document Mail Center – WO66-0609 10903 New Hampshire Avenue Silver Spring, MD 20993-0002

RE: Premarket Notification

To Whom It May Concern:

Enclosed in duplicate is the following information:

A. Purpose of Submission: <u>New Device</u>

B. Name and Address of the Third Party:



C. Name and Address of the Manufacturer:

Broncus Technologies, Inc. 1400 N Shoreline Blvd, Suite A8 Mountain View, CA 94043

D. Device Name

	LungPoint® Tools
Trade or Proprietary N	lame: (LungPoint Sheath, LungPoint Dilation Balloon)
Classification Name:	Bronchoscope and Accessories
Regulation Number:	21 CFR 874.4680
Recommendation:	Substantially Equivalent
Date Submission was	received by
(b)(4) Third Party Informati	on March 4, 2013


We have enclosed the following materials:

- E. Authorization Letter from the applicant (MAL-F-0006).
- F. Complete 510(k) application submitted by the applicant.
- G. Documented review of the 510(k) application (RPP-F-0012, RPP-F-14 and all correspondence and documents related to the review).
- H. Conflict of Interest Certification (RPP-F-0018)
- I. Certification (RPP-F-0020)

If you should have any questions regarding this submission please contact me at (b) (b)(4) Third or fax (b)(4) Third Party Information regarding this submission to (b)(4) Third Party Information


Sincerely. (b)(6)





Submission Certification

- 1. I certify that (b)(4) Third Party Information continues to meet the personnel qualifications and prevention of conflict of interest criteria reviewed by the FDA;
- In addition, I state that (b)(4) Third Party Information believes that statements made in the review are true and accurate to the best knowledge of (b)(4) Third Party Information
- 3. (b)(4) Third Party Information review is based on the 510(k) that is attached with the review; and
- (b)(4) Third Party Information understands that the submission of false information to the government is prohibited by 18 U.S.C. 1001 and 21 U.S.C. 33(q).





Conflict of Interest Declaration and Certification For the review of the 510(k) submission from

Applicant: Broncus Technologies, Inc.

Device Name or Model Name: LungPoint Tools

Initials
(b)(6) I have read and understand (b)(4) Third Party Information Confidentiality Procedure (COI-S-0023), regarding conflict of interests and the attachments accompanying the procedure and am aware of my responsibilities under them.
I have not been employed within the last twelve months by the firm who submitted the 510(k) for evaluation.
I did not charge fees contingent or based upon the recommendation for initial classification (SE decision).
I have not performed testing in connection with this specific device 510(k).
I understand that the Accredited Persons (AP) Program requires that the Accredited Person or any of its personnel involved in 510(k) reviews, which includes those who have authority over the review process, have no ownership or other financial interest in a device manufacturer or distributor that presents the appearance of a conflict of interest.
I do not participate in the design, manufacture or distribution of any medical device.
I do not provide consultative services to any device manufacturer or distributor regarding specific device 510(k) or participate in the preparation of 510(k.
I have not performed a 510(k) review where I have a personal relationship with the sponsor or the application correspondent.
(b)(6)
Signed:
Printed Nam

Date:

4 MAR 2613

(b)(4) Third Party Information

Questions ? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

Conflict of Interest Declaration and Certification For the review of the 510(k) submission from


Applicant: Broncus Medical, Inc.

LungPoint® Tools (LungPoint Sheath, Device Name or Model Name: LungPoint Dilation Balloon)

Initials I have read and understand (b)(4) Third Conflict of interest and Confidentiality Procedure (COI-S-0023), regarding conflict of interests and the attachments accompanying the procedure and am aware of my responsibilities under them. I have not been employed within the last twelve months by the firm who submitted the 510(k) for evaluation. I did not charge fees contingent or based upon the recommendation for initial classification (SE decision). I have not performed testing in connection with this specific device 510(k). I understand that the Accredited Persons (AP) Program requires that the Accredited Person or any of its personnel involved in 510(k) reviews, which includes those who have authority over the review process, have no ownership or other financial interest in a device manufacturer or distributor that presents the appearance of a conflict of interest. I do not participate in the design, manufacture or distribution of any medical device. I do not provide consultative services to any device manufacturer or distributor regarding specific device 510(k) or participate in the preparation of 510(k. I have not performed a 510(k) review where I have a personal relationship with the sponsor or the application correspondent.

	(b)(6)	
Signed:		
Printed Name:		
Date:	3 15 13	

Third Party Review Reviewer Memorandum



(b)(4) Revision (Effective August 1, 2011 Page 1 of 28

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118



I. Purpose and Submission Summary:

The 510(k) holder **Broncus Medical**, Inc. is submitting a Traditional Premarket Notification to introduce the **LungPoint® Tools (LungPoint Sheath, LungPoint Dilation Balloon)**, which are new medical devices, into interstate commerce. The Broncus Medical, Inc. LungPoint Software with FlexNeedle (previously cleared under K112051) is not included in the scope of this submission as the previously cleared LungPoint Software with FlexNeedle is not being modified in any way to allow for the use of these LungPoint Tools.

During the review of the submission dated March 1, 2013 one round of deficiencies was issued on March 20, 2013. Additional information dated April 3, 2013 was provided to respond to the deficiencies. A second round of deficiencies was issued on April 19, 2013. Additional information dated April 22, 2013 was provided to respond to the deficiencies. All deficiencies have been adequately addressed.

II. Administrative Requirements

	Yes	No	N/A
Indications for Use = Prescription Page Number: Attachment 2 of Additional Information dated April 3, 2013	x		
Truthful and Accuracy Statement Page Number: Page 19 of the original submission	x		
510(k) Summary Page Number: Attachment 2 of Additional Information dated April 22, 2013	x		
Standards Form (FDA Form 3654) Page Number: Appendix I of the original submission and Attachment 3 of Additional Information dated April 22, 2013	x		

III. Device Description

(4) Third Party Informa

	Yes	No	N/A
Is the device life-supporting or life sustaining?		Х	
Is the device an implant (implanted longer than 30 days)?		Х	
Does the device design use software?		Х	
Is the device sterile?	Х		
Is the device reusable (not reprocessed single use)?		Х	
Are "cleaning" instructions included for the end user?		Х	

The submission and additional information describes the LungPoint Sheath and LungPoint Dilation Balloon as follows:

Records processed under FOIA Request # 201



Accredited Person SE Documentation



Device Description

The updated device description is included in Attachment 6 of the Additional information dated April 3, 2013. The LungPoint Sheath and LungPoint Dilation Balloon are endoscopic tools, which are intended to be used as accessories to the Broncus Medical, Inc. LungPoint Software with FlexNeedle (previously cleared under K112051). The LungPoint Tools are described in more detail in the following subsections. The LungPoint Software (K112051) may be used with these tools or other commercially available tools. The LungPoint Tools have been specifically designed to be used with the LungPoint Software but can be used without the LungPoint Software, as well. The LungPoint Software with FlexNeedle (K112051) was not modified in any way to allow for the use of the LungPoint Tools and therefore is not included in the scope of this submission.

LunaPoint Sheath

The LungPoint Sheath is an endoscopic tool that is designed to be used with bronchoscopes to provide a working channel through which endoscopic tools, such as needles, dilation balloons, or other endoscopic devices may be introduced to the targeted area within the respiratory organs. The sheath enables physicians to easily access targets and allows for multiple approaches to the preselected target. It allows for the repeated placement of endoscopic tools to a specified lesion(s) during one procedure. The sheath is provided with (b)(4), which is

The sheath (b)(4) is inserted through a standard bronchoscope. It is then advanced to a target under the guidance of the LungPoint Software. Placement of the sheath may require the use of other endoscopic tools, (b)(4)

section. Once placed, other endoscopic tools

endoscopic devices, may be introduced through the sheath to the targeted area within the respiratory organs. (b)(4)

tools. The LungPoint Sheath is provided as a sterile device for single use (b)(4)

. A detailed drawing of the sheath is provided in Appendix B of the original submission (a portion of which is inserted below).

LungPoint Dilation Balloon

The LungPoint Dilation Balloon is an endoscopic tool that is used to dilate the target lung tissue of the bronchial tree. The dilation balloon is inserted through the bronchoscope or sheath and is used to dilate lung tissue. The balloon may be used (b)(4)

Revision (b) Effective August 1, 2011



Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

Records processed under FOIA Request # 20(b)(4) Third Placy Information

Accredited Person SE Documentation

(b)(4) as determined by the physician performing the bronchoscopy. The LungPoint Dilation Balloon (b)(4) , and is a sterile device for single use.

Balloon Size	Balloon	Maximum Catheter	Catheter Length/
(OD x length)	Pressure	Outer Diameter (OD)	working length
4mm x 6mm	10 atm	1mm	1430mm/975mm

When inflated at (b)(4) the balloon dimensions are (b)(4) . A detailed drawing of the LungPoint Dilation Balloon is provided in Appendix B of the original submission (a portion of which is inserted below).



Device Specifications

The following tables (extracted from Section 11 Device Description) outline the proposed **LungPoint Sheath and LungPoint Dilation Balloon** device performance characteristics and specifications.

Performance Specification	Objective	Applicable Standards	Acceptance Criteria (confidence/reliability parameters)
Dimensional Testing	To verify that the devices meet the product specification of the following Dimensional Tests: Inner diameter Outside diameter Sheath length Stylet tip diameter Stylet length	NA – all criteria are per the product specifications	(b)(4)
Joints/ Tensile Test	To verify that the joints of the sheath and stylet meet the tensile requirements of ISO 10555 and EN 1618	Method: EN1618:1997 Acceptance Criteria: adapted from ISO 10555- 1:1995	

LungPoint Sheath Specifications

b)(4) Third Party Information

Revision (b Effective A

Records processed under FOIA Request # 2015(b)(4) The Party Information

Accredited Person SE Documentation

Dec	Obia atian	A	A
Performance	Objective	Applicable	Acceptance Criteria
Specification		Standards	(b)(4)
Simulated Use	To verify that the sheath (b)(4) meet product specifications regarding use in a bronchoscope after repeated device insertions which constitute normal (b)(4)	NA – validated test fixture	(b)(4)
Radiopacity (b)(4)	NA	

LungPoint Dilation Balloon Specifications

Performance Specification	Objective	Applicable Standards	Acceptance Criteria (confidence/reliability parameters)
Dimensional Testing	To verify that the devices meet the product specification of the following Dimensional Tests: (b)(4)	NA – all criteria are per the product specifications	(b)(4)
Joints/ Tensile Test	To verify that the joints of the dilation balloon meet the tensile requirements of ISO 10555 and EN 1618	Method: EN1618:1997 Acceptance Criteria: adapted from ISO 10555- 1:1995	

b)(4) Third Party Information

Revision (b Effecti

4) I hird Party

Records processed under FOIA Request # 2015

b)(4) Third Party Information

Accredited Person SE Documentation

Performance Specification Simulated Use	Objective To verify that the balloon meets product specifications regarding	Applicable Standards NA – validated test fixture	Acceptance Criteria (confidence/reliability parameters) (b)(4)
	use in a bronchoscope after repeated device insertions which constitute normal use, in a test fixture (b)(4)		
Balloon burst pressure	To verify fatigue capabilities of the balloon and that balloon will not burst at or below maximum recommended burst pressure	NA	
Radiopacity	Verify that marker bands are visible under fluoroscopy (b)(4)	NA	

List of Materials

The following tables outline the lists of components and the patient contacting materials for the LungPoint Sheath and LungPoint Dilation Balloon.

LungPoint Sheath

Material Use	Material Description	Common Name
(b)(4)		
) Third Party Information		(b)(4) Third

Page 6 of 28 Questions ? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

Revision , Effective

b)(4) Manufacturing Information

LungPoint Dilation Balloon

(b)(4) Third Party Information

b)(4) Manufacturing Information

(b)(4) Third Party

Revision b Effecti

Page 7 of 28 Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

IV. Indications for Use

The proposed Indications for Use for the new device LungPoint Tools - LungPoint Sheath and LungPoint Dilation Balloon are as follows:

The LungPoint Tools are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs in a minimally invasive manner.

The Indications for Use listed above for the LungPoint Tools - LungPoint Sheath and LungPoint Dilation Balloon are the same as the Indications for Use stated in the Indications for Use form, the 510(k) Summary and the Instructions for Use.

V. <u>Predicate Device Comparison</u>

The updated Section 12 in Attachment 7 of the Additional Information dated April 3, 2013 and the responses in the Additional Information dated April 22, 2013 provide comparisons between the **LungPoint Sheath and LungPoint Dilation Balloon** submitted in this 510(k) and the predicate devices outlined below.

Predicate Device Trade Name	510(k) Submitter	510(k) Number	Reg # / Product Code
LungPoint Planning and Virtual Bronchoscopic Navigation (VBN) System (with FlexNeedle)	Broncus Technologies, Inc*	K112051, cleared to market on October 12, 2011 (and K093423, K091160 and K090095, by reference)	892.2050 / LLZ Picture Archiving and Communications System
inReach System	SuperDimension	K110093, cleared to market on February 11, 2011 (and K071473, K092365 and K102604, by reference)	892.1750 / JAK Computed Tomography X-Ray System
CRE Pulmonary Balloon Dilatation Catheter	Boston Scientific Corporation	K023337, cleared to market on November 18, 2002	874.4680 / KTI Bronchoscope and Accessories

*NOTE: Broncus Technologies, Inc changed their name to Broncus Medical, Inc in June 2012.

The LungPoint Tools (proposed regulation # 874.4680 and product code EOQ for Bronchoscope and Accessories), which are designed to be used as accessories to the LungPoint Software System that was cleared to market under K112051, are commonly used endoscopic tools. They are similar to other commercially available endoscopic tools but were designed to be used as accessories to the LungPoint Software.

b)(4) Third Party Information (b)(4) Third Party	
Revision <mark>(b)</mark> Effective August 1, 2011 Page 8 of 28	

DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

# :	$2^{(b)(4)}$ Third	Party Inform	ation	

Similarities and differences among these devices are outlined in the following comparison tables.

	i Logio in Roach System	LungBoint Software with FlowNeedle
	i Logic inReach System Components and Tools (K110093)	LungPoint Software with FlexNeedle (K112051) and LungPoint Tools (new device)
Intended Use (General	Enable the diagnosis and the treatment of lung cancer.	Enable the diagnosis and the treatment of lung cancer.
Purpose of Device – encompasses indication)	The inReach Tools and Components are designed to enable physicians to reach distant lung lesions in a	The LungPoint Tools are specifically designed for use with the LungPoint Software.
,	minimally invasive manner. They are specifically designed for use with the i.Logic System.	Endoscopic tools include FlexNeedle aspiration needle (previously cleared under K112051), sheath and dilation
	Endoscopic tools and components include biopsy forceps, cytology brushes, aspiration needles, locatable guide, guide catheter etc.	balloon (new devices).
Indications for Use for Software	The indications for use of the inReach System:	The LungPoint Tools are endoscopic tools used with bronchoscopes and
	Indicated for displaying images of the tracheobronchial tree to aid the physician in guiding endoscopic tools or catheters in the pulmonary tract and to enable marker placement within soft	intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs in a minimally invasive manner.
	lung tissue. It does not make a diagnosis and is not an endoscopic tool. Not for pediatric use.	The LungPoint Tools are being cleared for use with the LungPoint Software System. The indications for use of the LungPoint Software:
		Indicated for displaying images of the tracheobronchial tree to aid the physician in guiding endoscopic tools or catheters in the pulmonary tract and to enable marker placement within soft lung tissue. It does not make a diagnosis and is not an endoscopic tool. Not for pediatric use.
Radiopaque Markers	Yes	Yes
Delivery	Visual (via bronchoscope) and/or	Visual (via bronchoscope) and/or
Approach	Fluoroscopy	Fluoroscopy
Delivered	Flexible bronchoscope with minimum	Flexible bronchoscope with minimum
Through Single Lies	working channel of b mm	working channel of b mm
Single Use	Yes	Yes
Sterile	Yes	Yes

Comparison Table LungPoint Tools for use with LungPoint Software

b)(4) Third Party Information

Revision b Effective



Comparison Table LungPoint Sheath

	I-Logic inReach Extended Working Channel (K110093)	LungPoint Sheath (new device)
Indications for Use for Tools	There are no separate indications for use of the inReach System Tools and Components as they fall under the indications for the software: Indicated for displaying images of the tracheobronchial tree to aid the physician in guiding endoscopic tools or catheters in the pulmonary tract and to enable marker placement within soft lung tissue. It does not make a diagnosis and is not an endoscopic tool. Not for pediatric use.	The proposed indications for use for the LungPoint Tools are: The LungPoint Tools are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs in a minimally invasive manner.
Stylet Provided Stylet Use		
Catheter ID/OD	2.0/2.6 cm	2.0/2.6 cm
Catheter Length	975 mm	975 mm
Tensile Strength	Complies with ISO 10555 as it is the industry standard (assumed)	Complies with ISO 10555

Comparison Table LungPoint Dilation Balloon

	CRE Pulmonary Balloon	LungPoint Dilation Balloon
	Dilatation Catheter (K023337)	(new device)
Indications for	The cleared indications for use	The proposed indications for use for the
Use for Tools	for the CRE Pulmonary Balloon	LungPoint Tools are:
	Dilatation Catheter are: The CRE Pulmonary Balloon Dilatation Catheter is intended to	The LungPoint Tools are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint
	be used endoscopically to dilate strictures of the airway tree.	Software to aid in reaching a targeted area within the respiratory organs in a minimally invasive manner.
Balloon Size and	b)(4)	
Inflation (balloon		
pressure)		
Catheter Length		
Inflation Port		
Balloon Length		
Rated Burst		
Pressure		
Tensile Strength	Complies with ISO 10555 as it is	Complies with ISO 10555
	the industry standard (assumed)	

b)(4) Third Party Information

(b)(4) Third Party

Revision **(b)** Effective August 1, 2011 Page 10 of 28

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

Records processed under FOIA Request # 2^{(b)(4)} Third Party Information

Accredited Person SE Documentation

The following comparison information has been summarized from the submission and includes a comparison of the similarities and differences in the intended use/indications for use and the device specifications.

The LungPoint Sheath and LungPoint Dilation Balloon has the following similarities to the predicate devices:

- Same intended use (LungPoint Software with FlexNeedle plus LungPoint Sheath and LungPoint Dilation Balloon vs. inReach System Tools and Components)
- Same software indications for use (LungPoint Software with FlexNeedle vs. inReach System)
- Similar tools indications for use (LungPoint Sheath vs. inReach System Tools and Components)
- Similar tools indications for use (LungPoint Dilation Balloon vs. CRE Pulmonary Balloon Dilatation Catheter)
- Same target tissues (will reach distant lung lesions)
- Same operating principle (minimally invasive access through bronchoscopes)
- Similar fundamental scientific technology (sheaths and dilation balloon catheters utilize similar technologies, dimensions, and performance characteristics)
- Same dimensions for sheath/stylet tools (LungPoint Sheath vs. inReach System Tools and Components)
- · Same sterile condition and intended for single use only

The differences between the LungPoint Sheath and LungPoint Dilation Balloon and the predicate devices are as follows:

Several Design Differences

The development and design inputs of the proposed new device design were based on customer requirements in accordance with Broncus Medical's Design Control SOP. Both the new and predicate devices use the same method of action of dilating target tissue (b)(4) (b)(4) The key questions of safety and effectiveness relating to the dilation of target tissue in the lung are (b)(4)

(b)(4)

(b)(4) require more pressure to burst. The process of inflation/deflation of the LungPoint Dilation

Balloon remains the same.

(b)(4) Third Party Information

b)(4) Third Party

Revision **(b)** Effective August 1, 2011 Page 11 of 28

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118



(b)(4)

These differences in designs were supported by the full list of design verification testing conducted to demonstrate the new device complies with the Performance Specification Tables in the **Device Specifications** section on page 4 of this Review Memorandum. These test results including (b)(4) (b) were reviewed and found acceptable. Refer to **Section XI Performance Testing Bench** of this Review Memorandum for a summary of this design verification testing. These design differences have been adequately addressed.

Summary

Although there are some differences between the new devices and the predicate devices, these changes would not be expected to pose any additional or new risks beyond the predicate device designs since they performed equivalently in the biocompatibility testing and bench performance testing. The LungPoint Sheath and LungPoint Dilation Balloon met all the predetermined acceptance criteria of design verification as specified by applicable standards and FDA guidance. Therefore, these differences in the new device do not raise any new issues of safety and effectiveness, and the new LungPoint Sheath and LungPoint Dilation Balloon are equivalent in design and technological characteristics to the predicate devices.

Intended Use/Indications for Use

The proposed Indications for Use for the new device LungPoint Tools - LungPoint Sheath and LungPoint Dilation Balloon are as follows:

The LungPoint Tools are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs in a minimally invasive manner.

The proposed indications for use for the new device are very similar to the cleared indications

(b)(4)

Indicated for displaying images of the tracheobronchial tree to aid the physician in guiding endoscopic tools or catheters in the pulmonary tract and to enable marker placement within soft lung tissue. It does not make a diagnosis and is not an endoscopic tool. Not for pediatric use.

The proposed indications for use for the new device are very similar to the cleared indications for the CRE Pulmonary Balloon Dilatation Catheter predicate device (K023337) which are as follows:

The CRE Pulmonary Balloon Dilatation Catheter is intended to be used endoscopically to dilate strictures of the airway tree.

The LungPoint Sheath and LungPoint Dilation Balloon have the same indications for use as the predicate devices combined.

(b)(4) Third Party Information

b)(4) Third Party

Revision Effective August 1, 2011 Page 12 of 28

Questions: Contact - DAVODIANDCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118



Principle of Operation / Technological Characteristics


The LungPoint Sheath and LungPoint Dilation Balloon employ the same principles of operation and similar design and construction as the predicate devices.

Design / Materials / Specifications

The similarities and differences in physical designs and specifications between the new devices **LungPoint Sheath and LungPoint Dilation Balloon** and the predicate devices were outlined in the predicate device comparison tables above. Refer to **Section XI Performance Testing Bench** for a summary of the verification activities completed to demonstrate the design outputs of the new device met the design inputs and user need requirements. Since there are no significant differences in designs between the new devices and the predicate devices, there would not be any additional or new risks expected beyond the predicate device designs since they performed equivalently in the bench performance testing outlined below.

Performance

The new devices LungPoint Sheath and LungPoint Dilation Balloon submitted in this 510(k) and the predicate devices have very similar performance characteristics. To support the new device's design verification, the following tables outline the lists of performance testing that was completed using the new LungPoint Sheath and LungPoint Dilation Balloon. All LungPoint Sheath and LungPoint Dilation Balloon. All LungPoint Sheath and LungPoint Dilation Balloon. All LungPoint Sheath and LungPoint Dilation Balloon samples used for testing of each product attribute were representative of finished products. The new LungPoint Sheath and LungPoint Dilation Balloon met all of the design verification specifications/requirements as outlined in the Performance Specification Tables in the Device Specifications section on page 4 of this Review Memorandum. Section 18 and the test reports in Appendices H1 through H4 in the submission include summaries with detailed descriptions of the bench test methods, acceptance criteria, and conclusions demonstrating compliance with the product specifications. The passing results are summarized below. Refer to Section XI Performance Testing Bench for a more detailed summary of the verification activities completed.



each requirement in the following table has been modified from the FDA Guidance to rem b)(4)

b)(4) Third Party Information

)(4) Third Party

Page 13 of 28

Revision Effective August 1, 2011

Records processed under FOIA Request # 201


Accredited Person SE Documentation

	Requirements	Applicable? Yes/No	Rationale/Testing Conducted	
(b)	(4)			

(b)(4) Third Party Information

(b)(4) Third Party

Revision, Effective August 1, 2011 Page 14 of 28 Questions ? Contact PDAYCDKH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118



The following tests were performed on **LungPoint Sheath** devices that were subjected to the (b)(4), where applicable. All testing results met the predetermined acceptance criteria that were established in the test protocols. The following table summarizes the performance testing conducted to support the LungPoint Sheath performance. A detailed discussion is included in the test reports in Appendices H1 through H3 in the original submission.

Third Party Info

LungPoint Sheath Specifications and Performance Testing Results

Performance	Objective	Applicable	Acceptance Criteria	Results
Specification		Applicable Standards	Acceptance Criteria (confidence/reliability parameters)	
(b)(4)				Pass
				Pass
				Pass
				Pass
				l

f) find rang mornation	
	Revision Effective August 1, ^{ti} 2011
	Page 16 of 28
Questions? Contact FDA/CDRH/OCE	/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

(h)(1) T



The following tests were performed on LungPoint Dilation Balloon devices that were (b)(4) , where applicable. All testing results met the pre-determined acceptance criteria that were established in the test protocols. The following table summarizes the performance testing conducted to support the LungPoint Dilation Balloon performance. A detailed discussion is included in the test reports in Appendices H1 through H4 in the original submission.

LungPoint Dilation Balloon Specifications and Performance Testing Results

	Performance	Objective	Applicable	Acceptance Criteria	Results
	Performance Specification		Applicable Standards	Acceptance Criteria (confidence/reliability parameters)	
(b)(4	-)				Pass
					Pass
					rass
					Pass
					Pass
					rass

Revision Effect Page 17 of 28 ONS : CONTACT DAVED CHICOLD DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

b)(4)

The design verification performance test summary above demonstrates the LungPoint Sheath and LungPoint Dilation Balloon met all the predetermined product specifications and acceptance criteria as specified by the test protocols and/or customer inputs. Similar performance testing methods for intravascular catheter systems are outlined in the FDA Guidance entitled Guidance for Industry and FDA Staff - Non-Clinical Engineering Tests and Recommended Labeling for Intravascular Stents and Associated Delivery Systems and FDA recognized ISO 10555-1 standard. Therefore, the physical designs, materials, and specifications for the LungPoint Sheath and LungPoint Dilation Balloon have been tested by appropriate methods and no new issues of safety and effectiveness are raised based on the performance characteristics and testing of the new device.

Conclusion

The LungPoint Sheath and LungPoint Dilation Balloon have the same indications for use and the same intended use as the predicate devices. The LungPoint Sheath and LungPoint Dilation Balloon are equivalent in design and technological characteristics to the predicate devices. Based on the comparisons with the predicate devices above and the safety and performance testing results provided in the submission, the LungPoint Sheath and LungPoint Dilation Balloon do not raise new questions of safety or effectiveness.

VI. Labeling

The updated proposed Instructions for Use for the LungPoint Sheath and LungPoint Dilation Balloon device are provided in Attachment 3 of the Additional Information dated April 3, 2013 and Attachment 1 of the Additional Information dated April 22, 2013, respectively. The updated proposed pouch labels are provided in Attachment 8 of the Additional Information dated April 3, 2013. The updated LungPoint Software with Endoscopic Tools Instructions for Use Addendum DRAFT is provided in Attachment 9 of the Additional Information dated April 3, 2013.

The proposed labels include all the relevant warnings and cautions to the user as were in the predicate device labels such as Rx only, see packaging insert for explanation of symbols, for single use, sterilized using radiation, sterility guaranteed if package unopened and undamaged. (b)(4)

The Instructions for Use for this prescription only device includes all of the relevant warnings and precautions as the predicate device, such as, the prescription only caution statement, sterilization information including do not re-use or resterilize, for single patient use only statement, follow the bronchoscope manufacturer's Instructions for Use, the statement that this device should only be used by trained physicians, inspection of the device and package seal for any damage prior to use, the possible complications, (b)(4)

(b)(4) Third Party Information

Records processed under FOIA Request # 2015(b)(4) Third Party information.

Accredited Person SE Documentation

(b)(4)

The Instructions for Use includes the same indications for use as the Indications for Use form. The Device Description includes product configuration and dimensions table that outlines key information such as LungPoint Sheath (b)(4), LungPoint Dilation Balloon (b)(4)

s of the devices. The Contraindications section states there are no known contraindications. The Preparation for Use, Operating Instructions, and Storage sections provide detailed step by step directions to assist the user in proper use of the device along with any relevant warning or precautionary notes. The last note in the LungPoint Dilation Balloon Operating Instructions states, "NOTE: The

o)(4)

The Graphic Symbol Legend is included on the last page to explain all symbols used on the labels. No other specific claims are made which would raise questions of safety and effectiveness.

VII. Sterilization/Shelf Life/Reuse

The LungPoint Sheath and LungPoint Dilation Balloon device is provided sterile to the user and is labeled for single use only. The sterilization information included in Section 14 of the submission and the Additional Information dated April 3, 2013 is summarized in the table below.

	YES	NO
1. Sterilant:	Х	
a. Sterilization method description	Х	
(e.g., Steam, EtO, Radiation): Radiation (electron beam (E-beam))		
b. Dose, for radiation	X	
(e.g., 25 – 50 kGy): <mark>(b)(4)</mark>		
c. Sterilant residuals remaining on the device: N/A		Х
For EO, the maximum levels of residuals of EO and ethylene chlorhydrin that remain on		
the device (note: not to include ethylene glycol residual level because the recognized standard, "ANSI/AAMI/ISO 10993-7:1995 Biological Evaluation of Medical Devices –		
Part 7: Ethylene Oxide sterilization residuals," does not include measurement of		
ethylene glycol residuals);		
2. A description of the Validation Method for the sterilization cycle (not	X	
data): ISO 11137-1 :2006 and ISO 11137-2 :2012 using the method		
VDmax		
(Full citation of an FDA recognized standard is recommended		
(e.g., ANSI/AAMI/ISO 11135))		
3. Sterility assurance level (SAL): 10 ⁻⁶	X	
(e.g., 10 ⁻⁶ for all devices (except 10 ⁻³ for devices that contact intact		
skin))		
4. Is it labeled "Pyrogen Free"? No		Х
If so, a description of the method: N/A		Х
(e.g., LAL (<i>Limulus</i> Amebocyte Lysate test))		

b)(4) Third Party Information

(b)(4) Third Party

Revision () Effective August 1, 2011 Page 19 of 28 Records processed under FOIA Request # 201

Accredited Person SE Documentation

YES	NO
Х	
	YES X

Shelf Life

The **LungPoint Sheath and LungPoint Dilation Balloon** will have a shelf-life (b)(4). The data to support the shelf-life will be available prior to commercialization. The shelf life verification testing protocol includes (b)(4)

Shelf-life will be extended to up(b)(4) of representative samples (b)(4) (b)(4) (D)(4) (D)(4) (b)(4)

on the packaging. Testing of the LungPoint Sheath and LungPoint Dilation Balloon includes the following performance tests:

).

- Dimensional Testing (b)(4)
- Simulated Use Testing ((b)(4)
- Tensile (bb)(4) (b)(4) testing (b)(4)
- Balloon Fatigue and Burst Testing (b)(4)

The same performance tests will be conducted at each time interval as were tested for the design verification testing except for (b)(4) which is not affected by aging.

The submission includes all required details related to packaging, sterilization, and expiration dating. The device is intended for single use only.

VIII. Biocompatibility

The following tables outline the lists of components and the patient contacting materials for the **LungPoint Sheath and LungPoint Dilation Balloon**. For details on the colorants used in the **LungPoint Sheath**, refer back to the **List of Materials** Section starting on page 7 of this Review Memorandum. The colorant information provided by the sponsor was reviewed and found acceptable.

(b)(4) Third Party Information	Revision 6 Effective	(b)(4) Third Party Information
	E/DID at CDRH-FOISTATUS@fda.hhs.gov or 30	Page 20 of 28

Accredited Person SE Documentation

LungPoint Sheath

Material Use	Material Description	Common Name	Patient contact
(b)(4)			

LungPoint Dilation Balloon

	Material Use	Material Description	Trade Name	Patient contact
(b)	(4)			

Biocompatibility testing in accordance with FDA Blue Book Memo G95-1 and ISO 10993-1:2009 was required. The LungPoint Sheath and LungPoint Dilation Balloon are categorized as external communicating devices that will have limited exposure to mucosal membranes and lung tissues for <24 hr. Devices of this classification typically require the following tests: Cytotoxicity, Sensitization, Irritation/Intracutaneous Reactivity, and Acute Systemic Toxicity. The following tests were conducted: Cytotoxicity, Sensitization, Irritation/Intracutaneous Reactivity, and Acute Systemic Toxicity.

Biocompatibility studies were conducted by (b)(4) (b)(4) following ISO 10993 to ensure that the materials are safe to use for their intended duration and purpose. All parts coming into

b)(4) Third Party

Third Party Information	
	(b)(4) Third Pa
	Revision Effective August 1, 2011
	Page 21 of 28
	/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

Records processed under FOIA Request # 2015-(b)(4) Third Party Information

Accredited Person SE Documentation



direct or indirect contact with the patient (as outlined in the tables above) were tested. All testing was performed on finished, post-sterilized test articles. Test articles were extracted per the requirements in ISO 10993-12 (b)(4)

(b)(4)

All testing was completed in accordance with GLP regulations. Refer to the complete Biocompatibility Reports provided in Appendix G of the original submission. All reports were reviewed and found acceptable. Biocompatibility results summarized below and in these reports show that both devices meet all test requirements.

Test Name	Test Description	Results
Cytotoxicity	ISO 10993-5: In Vitro Cytotoxicity (L929 MEM Elution)	Pass (non-cytotoxic)
Sensitization	ISO 10993-10: Sensitization	Pass (non-sensitizing)
Intracutaneous Reactivity	ISO 10993-10: Irritation and Delayed-Type Hypersensitivity (Injection)	Pass (non-irritant)
Systemic Toxicity (Acute)	ISO 10993-11: Systemic Toxicity (Acute)	Pass (non-toxic)

The results of all biocompatibility testing were reviewed and found acceptable; no additional biocompatibility testing is required for the **LungPoint Sheath and LungPoint Dilation Balloon**. The LungPoint Sheath and LungPoint Dilation Balloon are considered safe and biocompatible for their intended use. No new questions related to safety and effectiveness are raised by the biocompatibility testing.

IX. Software

The LungPoint Tools (LungPoint Sheath and LungPoint Dilation Balloon) do not employ software. The Broncus Medical, Inc. LungPoint Software with FlexNeedle (previously cleared under K112051) is not included in the scope of this submission and the previously cleared LungPoint Software with FlexNeedle is not being modified to allow for the use of these LungPoint Tools. The existing risk analysis for the LungPoint Software was reviewed to assess whether guidance of the LungPoint Sheath and LungPoint Dilation Balloon adds any new hazards. No new risks were identified as the use of endoscopic tools, like the sheath and dilation balloon, is inherent to the design and intended use of the existing LungPoint Software. Additionally, no new verification and validation testing was performed as the LungPoint Software was not modified in any way to allow for the use of the LungPoint Tools.

X. Electromagnetic Compatibility and Electrical, Mechanical and Thermal Safety

The device is not electrically powered.

XI. Performance Testing – Bench

The new devices LungPoint Sheath and LungPoint Dilation Balloon submitted in this 510(k) and the predicate devices have very similar performance characteristics. To support the new device's design verification, the following tables outline the lists of performance testing that was completed using the new LungPoint Sheath and LungPoint Dilation Balloon. All LungPoint Sheath and LungPoint Dilation Balloon samples used for testing of each product attribute were

(b)(4) Third Party Information

απγ

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

Accredited Person SE Documentation

representative of finished products. The new LungPoint Sheath and LungPoint Dilation Balloon met all of the design verification specifications/requirements as outlined in the Performance Specification Tables in the **Device Specifications** section on page 4 of this Review Memorandum. Section 18 and the test reports in Appendices H1 through H4 in the submission include summaries with detailed descriptions of the bench test methods, acceptance criteria, and conclusions demonstrating compliance with the product specifications. The passing results are summarized below.

Refer to **Section V Predicate Device Comparison** for a more detailed summary of the design verification test methods that were utilized compared to the FDA Guidance entitled Guidance for Industry and FDA Staff - Non-Clinical Engineering Tests and Recommended Labeling for Intravascular Stents and Associated Delivery Systems (beginning on page 13 of this Review Memorandum).

The following tests were performed on **LungPoint Sheath** devices that were subjected to the required **b**(4), where applicable. All testing results met the predetermined acceptance criteria that were established in the test protocols. The following table summarizes the performance testing conducted to support the LungPoint Sheath performance. All samples sizes for this table were **(D**)(4). A detailed discussion is included in the test reports in Appendices H1 through H3 in the original submission.

Performance Specification	Objective	Applicable Standards	Acceptance Criteria (confidence/reliability parameters)	Results
Dimensional Testing	To verify that the devices meet the product specification of the following Dimensional Tests: Inner diameter Outside diameter 0(4)	NA – all criteria are per the product specifications	(b)(4)	
Joints/ Tensile Test	To verify that the joints of the sheath and stylet meet the tensile requirements of ISO 10555 and EN 1618	Method: EN1618:1997 Acceptance Criteria: adapted from ISO 10555- 1:1995		

LungPoint Sheath Design Verification Testing Summary



Records pro	cessed under FOIA Req	uest # 2015 <mark>(b)(4</mark>)	Third Party Information	
Accredited P	erson			
SE Documen				
Performance Specification	Objective	Applicable Standards	Acceptance Criteria (confidence/reliability parameters)	Results
Simulated Use	(b)(4)	NA – validated test fixture	parameters) (b)(4)	
Radiopacity	(b)(4)	NA		

The following tests were performed on LungPoint Dilation Balloon devices that were subjected to the required (b)(4), where applicable. All testing results met the pre-determined acceptance criteria that were established in the test protocols. The following table summarizes the performance testing conducted to support the LungPoint Dilation Balloon performance. All samples sizes for this table (b)(4) A detailed discussion is included in the test reports in Appendices H1 through H4 in the original submission.

Accredited Person SE Documentation

Performance Specification	Objective	Applicable Standards	Acceptance Criteria (confidence/reliability parameters) (b)(4)	Results
Dimensional Testing	(b)(4)	NA – all criteria are per the product specifications	(b)(4)	
Joints/ Tensile Test		Method: EN1618:1997 Acceptance Criteria: adapted from ISO 10555- 1:1995		
Simulated Use		NA – validated test fixture		

LungPoint Dilation Balloon Design Verification Testing Summary

Records processed under FOIA Request # 2015

Accredited Person SE Documentation

Performance Specification	Objective	Applicable Standards	Acceptance Criteria (confidence/reliability parameters)	Results
Balloon burst pressure	(b)(4)	NA	-(b)(4)	
Radiopacity	(b)(4) (b)(4) n.	NA		

The design verification performance test summary above demonstrates the LungPoint Sheath and LungPoint Dilation Balloon met all the predetermined product specifications and acceptance criteria as specified by the test protocols and/or customer inputs. Similar performance testing methods for intravascular catheter systems are outlined in the FDA Guidance entitled Guidance for Industry and FDA Staff - Non-Clinical Engineering Tests and Recommended Labeling for Intravascular Stents and Associated Delivery Systems and FDA recognized ISO 10555-1 standard. Therefore, the physical designs, materials, and specifications for the LungPoint Sheath and LungPoint Dilation Balloon have been tested by appropriate methods and no new issues of safety and effectiveness are raised based on the performance characteristics and testing of the new device.

Conclusion: The bench test results support the conclusion that the performance of the **LungPoint Sheath and LungPoint Dilation Balloon** is similar to the performance of the predicate devices. No new questions of safety and effectiveness are raised.

XII. Performance Testing – Animal

This submission does not include animal testing or data.

XIII. Performance Testing – Clinical

This submission does not include human clinical testing or data.

(b)(4) Third Party Information (b)(4) Third Party Information Revision Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118 Accredited Person SE Documentation

XIV. Substantial Equivalence Discussion

	Yes	No	
1. Same Indication Statement?	Х		If YES = Go To 3
2. Do Differences Alter The Effect Or Raise New Issues of Safety Or Effectiveness?			If YES = Stop NSE
3. Same Technological Characteristics?	Х		If YES = Go To 5
4. Could The New Characteristics Affect Safety Or Effectiveness?			If YES = Go To 6
5. Descriptive Characteristics Precise Enough?		Х	If NO = Go To 8 If YES = Stop SE
6. New Types Of Safety Or Effectiveness Questions?			If YES = Stop NSE
7. Accepted Scientific Methods Exist?			If NO = Stop NSE
8. Performance Data Available?	Х		If NO = Request Data
9. Data Demonstrate Equivalence?	Х		Final Decision: SE

Note: Document the decision path by marking the arrows followed on the FDA flowchart.

Please complete the following table and answer the corresponding questions. "Yes" responses to questions 2, 4, 6, and 9, and every "no" response requires an explanation.

- 1. Explain how the new indication differs from the predicate device's indication:
- 2. Explain why there is or is not a new effect or safety or effectiveness issue:
- 3. Describe the new technological characteristics:
- 4. Explain how new characteristics could or could not affect safety or effectiveness:
- 5. Explain how descriptive characteristics are not precise enough:

The submission includes the descriptive characteristics but the performance testing is needed to support substantial equivalence and demonstrate the similarities between the new devices and the predicate devices.

- 6. Explain new types of safety or effectiveness question(s) raised or why the question(s) are not new:
- 7. Explain why existing scientific methods cannot be used:
- 8. Explain what performance data is needed:

b)(4) Third Party Information

Revision (b) Effective Argust 1, 2011

Page 27 of 28

Records processed under FOIA Request # 201(b)(4) Third Party Information

Accredited Person SE Documentation



9. Explain how the performance data demonstrates that the device is or is not substantially equivalent:

As the reviewer of this submission I have reviewed the instructions for use, the submitter's description of the devices and compared this information against the information concerning the predicate devices that was provided by the submitter. The specifications for the predicate devices and the new devices have been compared. They are very similar. The predicate device comparison tables demonstrate the similarities and differences between the new devices and predicate devices. The submission includes performance and biocompatibility testing which demonstrates the new devices and the predicate devices possess similar performance characteristics. The labeling included in the submission was reviewed and found to be very similar to the predicate labeling. There are no new questions of safety and effectiveness raised during this review.

Based upon the above summary, a substantially equivalent decision is recommended.

XV. Deficiencies

During the review of the submission dated March 1, 2013 one round of deficiencies was issued on March 20, 2013. Additional information dated April 3, 2013 was provided to respond to the deficiencies. A second round of deficiencies was issued on April 19, 2013. Additional information dated April 22, 2013 was provided to respond to the deficiencies. All deficiencies have been adequately addressed.

XVI. Contact History

All correspondence is included in the submission.


XVII. Recommendation

Regulation Number:	21 CFR 874.4680
Regulation Name:	Bronchoscope and Accessories
Regulatory Class:	II
Product Code:	EOQ

(b)(4) Third Party Information

Revision b Effective A

Records processed under FOBr Brougst # 2015 78 fog Res 1996 by CDRH on 02-01-2016 LungPoint Tools (LungPoint Sheath, LungPoint Dilation Balloon)



* 510(k) Submissions compare new devices to marketed devices. FDA requests additional information if the relationship between marketed and "predicate" (pre-Amendments or reclassified post-Amendments) devices is unclear.

** This decision is normally based on descriptive information alone, but limited testing information is sometimes required.

*** Data maybe in the 510(k), other 510(k)s, the Center's classification files, or the literature.

This page represents the additional information provided by the sponsor.

Please see the enclosed Sponsor eCopy in Folder labeled as VOL_003_Response to 19Apr13.

MARK	JOB
------	-----

From:	(b)(6) (b)(4) Third Party
Sent:	Friday, April 19, 2013 2:17 PM
To:	(b)(4) Third Party
Cc:	(b)(6)
Subject:	Review of Responses dated April 3, 2013
Attachments:	Broncus Medical LungPoint Tools Deficiencies 19 Apr 2013.doc

Dear (b)(6)

We have completed our review of the Responses dated April 3, 2013 to (b)(4) Third Party Information

(b)(4) These deficiencies are described in the attached Record of Deficiencies from the Substantive Review. Your submission has been placed on hold until we receive your response to the deficiencies.

If you have any questions, please do	not hesitate to contact me.	You may reach me at	(b)(4)	or this email
address: (b)(4)		-		

Sincerely,

(b)(6)

Reviewer

(b)(4)

Record of Deficiencies From Substantive Review

Device Name or Model Name: LungPoint® Tools (LungPoint Sheath, LungPoint Dilation Balloon)

Date: April 19, 2013

Please provide the following items as required according to the 510(k) checklist, the FDA Guidance for Industry and FDA Staff on Format for Traditional and Abbreviated 510(k)s issued August 12, 2005, and the Draft Guidance on Submission and Review of Sterility Information in Premarket Notification (510(k)) Submissions for Devices Labeled as Sterile issued December 12, 2008.

We request that you provide detailed item-by-item responses to this request for additional information. The suggested format is to restate the question and/or comment, state your response to the request/recommendation, and reference the location of the appropriate supporting documentation in the response documents.

.

(b)(4) Revision BEffective May 30, 2003 DRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118 This page represents the additional information provided by the sponsor.

Please see the enclosed Sponsor eCopy in Folder labeled as VOL_003_Response to 20Mar13.

MARK JOB

(b)(6) (b)(4)
Wednesday, March 20, 2013 12:08 PM
(b)(4)
(b)(6)
Questions on Broncus Medical LungPoint Tools
Broncus Medical LungPoint Tools Deficiencies 20 Mar 2013.doc

Dear (b)(6)

We have completed both the administrative and the substantive review according to the 510(k) checklist, the FDA Guidance for Industry and FDA Staff on Format for Traditional and Abbreviated 510(k)s issued August 12, 2005, and the Draft Guidance on Submission and Review of Sterility Information in Premarket Notification (510(k)) Submissions for Devices Labeled as Sterile issued December 12, 2008.

(b)(4) These deficiencies are described in the attached Record of Deficiencies from the Substantive Review. Your submission has been placed on hold until we receive your response to the deficiencies.

If you have any questions, please do not hesitate to contact me. You may reach me at (b)(4) or this email address: (b)(4)

Sincerely,

<mark>(b)(6)</mark> Reviewer

(b)(4)

Record of Deficiencies From Substantive Review

Device Name or Model Name: LungPoint® Tools (LungPoint Sheath, LungPoint Dilation Balloon)

Date: March 20, 2013

Please provide the following items as required according to the 510(k) checklist, the FDA Guidance for Industry and FDA Staff on Format for Traditional and Abbreviated 510(k)s issued August 12, 2005, and the Draft Guidance on Submission and Review of Sterility Information in Premarket Notification (510(k)) Submissions for Devices Labeled as Sterile issued December 12, 2008.

We request that you provide detailed item-by-item responses to this request for additional information. The suggested format is to restate the question and/or comment, state your response to the request/recommendation, and reference the location of the appropriate supporting documentation in the response documents.

(b)(4)

Date: March 4, 2013

(b)(6)

Broncus Technologies, Inc. 1400 N Shoreline Blvd, Suite A8 Mountain View, CA 94043

Re: LungPoint Tools

Dear (b)(6)

This letter is to acknowledge on March 4, 2013, (b)(4) (b)(4) received the 510(k) dated March 1, 2013 for the LungPoint Tools.

We will keep you informed as the review progresses.

If you have any questions, please do not hesitate to contact me. You may reach me (b)(4) or email at (b)(4)

Sincerely,



Part Acceptance / Non-acceptance

1. Accredited Person:

Name:	(b)(4) Third Party Information	٦		 -
Address				 -
Contact:	(b)(6)			
Telephone	(b)(4) Third Party	Fax:	(b)(4) Third P t	

2. Foreign Accredited Person, Specify a Domestic Correspondent:


Name:	N/A		
Address			
Contact:			
Telephone:		Fax:	

3. 510(k) Owner (Applicant, Manufacturer, other persons preparing 510(k))

Name:	Broncus Medical,	Inc.				
Address	1400 N. Shoreline Boulevard, Suite A8					
	Mountain View, CA 94043					
Contact:	(b)(6) , R	egulatory Affairs Consultant				
Telephone:	650-428-1600	Fax: 650-428-1542 (b)(4)				

STOP!

Before completing items 4 to 9 below, complete pages 3-6 of this document.



4. Device Name:

LungPoint® Tools (LungPoint Sheath, LungPoint Dilation Trade or Proprietary Name: Balloon) Classification Name: Bronchoscope and Accessories; product code EOQ

- 5. CFR Classification Citation: 21 CFR 874.4680 (see 21 CFR 862 through 892)
- 6. Classification Panel: Ear, Nose, and Throat
- 7. Based on my completion of this document, I recommend that this 510(k):



- Be accepted for substantive review and I have notified the 510(k) owner using RPP-F-0016.
- Π


Not be accepted for substantive review and I have listed the deficiencies on RPP-F-0016.

8. Primary Reviewer



3/15/13

Date



<u>4/29/13</u> Date



Print Name



CI	necklist Questions:	YES	NO	Instructions
1.	a). Is the device one that FDA has determined as being acceptable for third party review?			If NO, telephone DSMA for instructions. STOP REVIEW
1	b). Have you confirmed that the manufacturer has not engaged in forum shopping?			If NO, telephone DSMA for instructions. STOP REVIEW
2.	Is the device trade or proprietary name included?			If NO, note deficiency on RPP-F-0013.
3.	Is the device common or usual name included?			If NO, note deficiency on RPP-F-0013.
4.	Is the device classification name, class of the device, and regulation number (21 CFR 874.4680; product code EOQ) included?			If NO, note deficiency on RPP-F-0013.
5.	Is the classification panel included?			If NO, note deficiency on RPP-F-0013.
6.	Has the applicant complied with Section 514 of the Act? (Section 514 relates to performance standards for class II devices. At this time, there are no 514 standards. Therefore, your answer should be yes.)	\boxtimes		lf NO, note deficiency on RPP-F-0013.
7.	Does the submission include proposed labels, labeling, and advertisements (if available) that describe the device, its intended use, and directions for use (ODE Guidance Memorandum #G91-1)?			lf NO, note deficiency on RPP-F-0013.
8.	Does the submission contain the "Indications for Use" form?			If YES, indicate page number <u>17</u> . If NO, note deficiency on RPP-F-0013.

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016 Acceptance Checklist

2

Checklist Questions:	YES	NO	Instructions
9. Does the submission contain an acceptable <u>510(k) Summary</u> of Safety and Effectiveness (per 21 CFR 807.92) OR an acceptable <u>510(k) Statement</u> (per 21 CFR 807.93) that safety and effectiveness information will be made available to any person upon request?			If YES, indicate page number <u>Appendix A</u> . If NO, note deficiency on RPP-F-0013.
10. Does the submission contain photographs of the device if applicable?			If NO, note deficiency on RPP-F-0013.
11. Does the submission contain drawings for the device with dimensions and tolerances if applicable?			If NO, note deficiency on RPP-F-0013.
12. Does the submission identify the device to which equivalence is claimed?			
13. If the answer to question 12 is YES, did the applicant identify:			
a. Predicate device (referred to as marketed device)?		\boxtimes	
b. Legally marketed device (referred to as marketed device)?	\boxtimes		Note deficiency on RPP-F-0013.
Note: A predicate device is a device that was legally in commercial distribution in the U.S. on or before May 28, 1976 (referred to as a pre- amendments device) or a device that was marketed after May 28, 1976 (referred to as a post amendments device) that was reclassified from class III to class I or II. A marketed device can be a predicate device but is most often a device that FDA has determined is SE to another marketed device (21 CFR 807.92(a)3). IT IS YOUR RESPONSIBILITY TO MAKE SURE THAT THE PREDICATE DEVICE OR LEGALLY MARKETED DEVICE IDENTIFIED IS LEGITIMATE. If it is not, the review must STOP. Telephone DSMA for assistance.			List all 510(k) control numbers: <u>K112051 (SE 10/12/2011)</u> <u>K110093 (SE 02/11/2011)</u> <u>K023337 (SE 11/18/2002)</u>

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016 Acceptance Checklist

.

Checklist Questions:	YES	NO	Instructions	
14. Does the submission contain information about the marketed device(s) identified in questions 12 and 13 above to which equivalence is claimed, including labeling and a description of the device?			If NO, note deficiency on RPP-F-0013.	
15. Does the submission contain a statement/comparison of similarities and/or differences between the new device and the marketed device? (The new device that is the subject of this 510(k) can be either a new device or a modification to the existing device.)			If NO, note deficiency on RPP-F-0013.	
16. Does the submission contain the Truthful and Accurate Statement (per 21 CFR 807.87(j))?			If YES, indicate page number <u>19</u> . If NO, note deficiency on RPP-F-0013.	
17. Does the submission contain the submitter's name, address, contact person, telephone number, and fax number?	\boxtimes		If NO, note deficiency on RPP-F-0013.	
18. If there is a representative or consultant, does the submission contain their name, address, contact person, telephone number, and fax number?			If NO, note deficiency on RPP-F-0013.	
19. Does the submission contain a table of contents with pagination?	\boxtimes		If NO, note deficiency on RPP-F-0013.	
20. If the submitter has a manufacturing facility (contract or owned), and/or a sterilization facility (contract or owned), is the address(es) contained in the submission?			If NO, note deficiency on RPP-F-0013.	
21. Does the submission contain a comparison table of the new device to the marketed device?			If NO, note deficiency on RPP-F-0013.	
22. Does the submission contain information about the action taken to comply with voluntary standards?			If NO, note deficiency on RPP-F-0013.	
(b)(4) Revision(b, Effective 01 October 03) Page 5 of 6 FOISTATUS@fda.hhs.gov or 301-796-8118				

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016 Acceptance Checklist (b)(4)

Checklist Questions:	YES	NO	Instructions
23. Does the submission contain performance data (can be bench or animal but not clinical), i.e.:			
Is there performance data for the marketed device? a. Bench testing? b. Animal testing?			If NO, note deficiency on RPP-F-0013. <u>Predicate</u> device bench and animal testing is not required.
Is there performance data for the new device? a. Bench testing? b. Animal testing?			If NO, note deficiency on RPP-F-0013. <u>New device</u> animal testing is not required.
24. If the device is labeled as sterile, does the submission contain sterilization data?			If NO, note deficiency on RPP-F-0013.
25. Does the device incorporate a computer or computer software?		\boxtimes	
a. If YES, is there information about the hardware?			If NO, note deficiency on RPP-F-0013. <u>No software</u>
b. If YES, is there information about the software?			
26. a) Is there a specific guidance document for this type of device? Title:			If YES, continue review with checklist from the specific guidance document and return to question 27.
			If NO, proceed to question 26 b).
26 b) Contact the appropriate ODE Branch Chief to obtain information for reviewing this type of device. Has a summary of this discussion been documented?			If YES, answer question 27. If NO, do not proceed to question 27; stop review until summary completed.
27 Is this 510(k) sufficiently complete to allow substantive review?			If YES, continue review using specific guidance document or if no specific guidance document, continue the review using documentation forms. If NO, note deficiency on RPP-F-0013.

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016



Broncus Medical, Inc. 1400 N. Shoreline Blvd. Suite A8 Mountain View, CA 94043 www.broncus.com (650) 428-1600 (650) 428-1542 fax

March 1, 2013



Re: Authorization for Accredited Person Review of 510(k)

Dear (b)(6)

Enclosed is the Traditional 510(k) for Broncus' LungPoint Tools.

We at Broncus Medical, Inc hereby authorize (b)(4) to submit the enclosed 510(k) to the Food and Drug Administration (FDA) on our behalf, discuss its contents with the FDA, and function as the Accredited Person to perform the third party review.

We certify that we have not established a contract with another Accredited Person to perform the review of this 510(k) submission.

We accept the quote for the Traditional 510(k) review services, including the (b)(4) (b)(4) Terms and Conditions.

Additionally, a completed Refuse to Accept Checklist is attached for your reference. Should you have any questions, please contact me at (b)(4) (b)(4) or by e-mail at (b)(4)

Sincerely,

b)(6)

Consulting Regulatory Affairs Director

Enclosure

TRADITIONAL 510(k)

LungPoint® Tools

Broncus Medical, Inc. 1400 N. Shoreline Boulevard Suite A8 Mountain View, CA 94043

Table of Contents

1. MDUFMA COVER SHEET	6
2. CDRH PREMARKET REVIEW SUBMISSION COVER SHEET	7
3. 510(K) COVER LETTER	14
4. INDICATIONS FOR USE STATEMENT	16
5. 510(K) SUMMARY	18
6. TRUTHFUL AND ACCURACY STATEMENT	19
7. CLASS III SUMMARY AND CERTIFICATION	20
8. FINANCIAL CERTIFICATION OR DISCLOSURE STATEMENT	21
9. DECLARATIONS OF CONFORMITY AND SUMMARY REPORTS	22
10. EXECUTIVE SUMMARY	22
10.1. Background and Rationale	23
10.2. Description of Device	23
10.3. Device Comparison	23
10.4. Summary of Performance Testing	23
11. DEVICE DESCRIPTION	24
11.1. LungPoint Sheath	24
11.1.1. Materials	24
11.1.2. Performance Specification	26
11.2. LungPoint Dilation Balloon	28
12. SUBSTANTIAL EQUIVALENCE DISCUSSION	31
12.1. Predicate Devices	31
12.2. Comparison to the Predicate Devices	31
12.3. Substantial Equivalence Discussion	33
13. PROPOSED LABELING	36
14. STERILIZATION/SHELF LIFE	37
14.1. Overview	37
14.2. Sterilization Process	38
14.3. Dose Audit	38
14.4. Sterilization Validation and Bioburden	38

14.5. Dose Mapping Study for Product Release	39
14.6. Pyrogenicity	39
14.7. Shelf-Life	39
15. BIOCOMPATIBILITY	40
15.1. Overview	40
15.2. Materials in Contact with Patients	40
15.3. Biocompatibility Summary	41
16. SOFTWARE	43
17. ELECTROMAGNETIC COMPATIBILITY/ELECTRICAL SAFETY	44
18. PERFORMANCE TESTING – BENCH	45
18.1. Overview	45
18.2. LungPoint Dilation Balloon	48
18.3. LungPoint Sheath	50
19. PERFORMANCE TESTING – ANIMAL	54
20. PERFORMANCE TESTING - CLINICAL	55

List of in-text figures

	Figure 1: 510(k) "Substantial Equivalence" Decision-Making Flowchart	.35
b)(4)		
		.50
		.50

List of in-text tables

(b)(4)	22
	24
	26
	28
	29
	32
	33
	40
	41
	41
	42
	48
	52

Appendices

Appendix A: 510(k) Summary Appendix B: Device Drawings Appendix C1: Supporting Information for (b)(4) Appendix C2: Supporting Information for (b)(4) Appendix D: Labeling Appendix E: Predicate 510(k) Summaries Appendix F: Shelf-life Protocol Appendix G: Biocompatibility Repor(b)(4) Appendix H1 (b)(4) Appendix H2 Appendix H3 Appendix H4 Appendix H5: Radiopacity Summary Appendix I: FDA Forms 3654

Abbreviation or Terminology	Definition
FMECA	Failure modes and effects criticality analysis
HID	Hazard identification
QA	Quality assurance
R&D	Research and development
V&V	Verification and validation

ABBREVIATIONS & TERMINOLOGY

(The remainder of this page is intentionally left blank.)

1. MDUFMA COVER SHEET

This section is not applicable since the submission is being submitted through (b)(4) b)(4), a third party reviewer, who is listed on the FDA website (as of 2/28/2013) as an accredited person for premarket submission review.

(The remainder of this page is intentionally left blank.)

2. CDRH PREMARKET REVIEW SUBMISSION COVER SHEET

FDA Form 3514 is included on the following page.

(The remainder of this page is intentionally left blank.)

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016												
DEPARTMENT OF HEALTH AND HUMAN SERVICES Form Approval FOOD AND DRUG ADMINISTRATION OMB No. 0910-0120												
CDRH PREMARKET REVIEW SUBMISSION COVER SHEET Expiration Date: December 31, 2013 See PRA Statement on page 5.												
Date of Submission	User Fee Payment	ID Number		FDA Submiss	ion Docume	ent Numbe	er (if known)					
03/01/2013	NA											
SECTION A TYPE OF SUBMISSION												
PMA Original Submission Premarket Report Modular Submission Amendment Report Report Report Amendment Licensing Agreement	PMA & HDE Supplement Regular (180 day) Special Panel Track (PMA Only) 30-day Supplement 30-day Notice 135-day Supplement Real-time Review	PD Original PI Notice of C	DP Completion	510(k)	l (Complete age 5)	Request for Feedback Pre-Submission Informational Meeting Submission Issue Meeting Day 100 Meeting Agreement Meeting Determination Meeting Study Risk Determination						
	Amendment to PMA & HDE Supplement Other		ation Datition	Evolución of Automotio		Othe	er (specify):					
IDE	Humanitarian Device Exemption (HDE)	Class II Exem		Evaluation of Automatic Class III Designation (De Novo) Original Submission Additional Information		Other Submission 513(g) Other (describe submission):						
Have you used or cited Stan		Yes No		please complete Se	ection I, Pag	e 5)						
SECTION B Company / Institution Name	SUBM	IITTER, APPLI		ONSOR Registration Number	(if known)							
Broncus Medical, Inc			3007867778	i ogioti di controli i ogio	(1110111)							
Division Name (if applicable)			Phone Number (including area code)									
			650-428-1600									
Street Address			FAX Number (including area code)									
1400 N Shoreline Ave Suite A	8		650-428-1542									
City			State / Provinc	e	ZIP/Postal	I Code	Country					
Mountain View			CA		94043		USA					
Contact Name		,					1					
(b)(6)												
Contact Title	· · · · · · · · · · · · · · · · ·	·	Contact E-mail	Address								
Regulatory Affairs Consultant			(h)(4)									
SECTION C Company / Institution Name	APPLICATION CORRES	SPONDENT (e.	g., consultar	it, if different fror	n above)							
Division Name (if applicable)			Phone Number (including area code)									
Street Address	FAX Number (including area code)											
City	State / Province ZIP Code		Country									
Contact Name												
Contact Title	Contact E-mail Address											
FORM FDA 3514 (1/13) Page 1 of 5 Pages												

PSC Publishing Services (301) 443-6740 EF

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301P23968818f 55

Records	processed under	FOIA Request #	ŧ 2015-7253; I	Released by	CDRH on	02-01-2016
---------	-----------------	----------------	----------------	-------------	---------	------------

SECTION D1 REA	ASON FOR APPLICATION - PMA, PDP, OR H	HDE
 New Device Withdrawal Additional or Expanded Indications Request for Extension Post-approval Study Protocol Request for Applicant Hold Request for Removal of Applicant Hold Request to Remove or Add Manufacturing Site Process change: Manufacturing Packaging Sterilization Other (specify below) Response to FDA correspondence: 	Change in design, component, or specification: Software / Hardware Color Additive Material Specifications Other (specify below) Labeling change: Indications Performance Characteristics Shelf Life Trade Name Other (specify below)	Location change: Manufacturer Sterilizer Packager Report Submission: Annual or Periodic Post-approval Study Adverse Reaction Device Defect Amendment Change in Ownership Change in Correspondent Change of Applicant Address
Other Reason (<i>specify</i>):		<u>I</u>
SECTION D2	REASON FOR APPLICATION - IDE Change in: Correspondent / Applicant Design / Device Informed Consent Manufacturer Manufacturing Process Protocol - Feasibility Protocol - Other Sponsor Report submission: Current Investigator Annual Progress Report Site Waiver Report Final	 Response to FDA Letter Concerning: Conditional Approval Deemed Approved Deficient Final Report Deficient Progress Report Deficient Investigator Report Disapproval Request Extension of Time to Respond to FDA Request Hearing
Other Reason (specify):	L	
SECTION D3	REASON FOR SUBMISSION - 510(k)	.4
New Device	Additional or Expanded Indications	Change in Technology
Other Reason (specify): FORM FDA 3514 (1/13)		Page 2 of 5 Pages

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301P2696916f 55

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016

SE	CTION E		ADDITION	NAL INFORMATION	N ON 5	10(<u>K) SU</u>	BMIS	SIO	NS	····	
Pro	oduct codes of devices	s to which substantial		is claimed						Summary of, or safety and effect	statement concerning, stiveness information
1	EOQ	2		3	4 🛛 🔀 510 (k) summary a) summary attached	
5		6		7 8				510 (k) statement		
	nformation on devices to which substantial equivalence is claimed (if known)										
	510	D(k) Number		Trade or Proprie	tary or I	Model Name				Man	ufacturer
1	K112051		1	LungPoint Software			1	Broncus Medical, Inc			
2	K110093		2	inReach System 2 SuperDimension (Covide					idien)		
3	K023337		3	CRE Pulmonary Balloc	on Dilata	ition Cathete	r	3	Bost	ton Scientific Corp	oration
4			4					4			
5			5					5			
6			6					6			
	CTION F		OUCT INFO	DRMATION - APPL	CATIO	ON TO AL	L AP	PLIC	ATI	ONS	
		or classification name									
Bı	onchoscope and Acce	ssories									
r	Trade or Proprietary	or Model Name for Th	is Device				М	odel N	lumb	er	
1	LungPoint Tools (Lu	IngPoint Sheath and Lu	ngPoint Dila	ation Balloon)			1				····
2							2				
3							3				
4							4				
5							5				
FD.	A document numbers	of all prior related sub	missions (reg	garaless of outcome)	4			5			6
7		8	9		10			11			12
Dat	ta Included in Submis		ratory Testing	a 🗆	nimal T	rials				Human Trials	L
SF	CTION G						LLA	PPL			
Pro	oduct Code	C.F.R. Section (if app					ce Cla				
	DQ	21 CFR 874.4680					Class	1	\boxtimes	Class II	د •
	Classification Panel Bronchoscope and Accessories Class III Unclassified										
Tł	ications (from labeling ne LungPoint Tools are ithin the respiratory or	e endoscopic tools used	l with bronche	oscopes and intended to	be used	as accessori	es to th	ne Lur	ıgPoir	nt Software to aid i	n reaching a targeted area
FO	RM FDA 3514 (1/	13)									Page 3 of 5 Pages

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 3(Page-&016f 55

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016

<i>Note:</i> Submission of the in need to submit device est	nformation entered in Section H d ablishment registration.	oes not affect the	FDA Document Number (if kno	(האס
SECTION H	MANUFACTURING /	PACKAGING / ST	FERILIZATION SITES REL	ATING TO A SUBMISSION
Original	Facility Establishment Identifier	(FEI) Number	Manufacturer	Contract Sterilizer
Add Delete			Contract Manufacturer	Repackager / Relabeler
Company / Institution Nan	he		Establishment Registration Nu	mber
Division Name (if applicat	le)		Phone Number (including area	a code)
Street Address			FAX Number (including area c	ode)
City			State / Drevin as	
City			State / Province	ZIP Code Country
Contact Name		Contact Title		Contact E-mail Address
Original	Facility Establishment Identifier ((FEI) Number	Manufacturer	Contract Sterilizer
Add Delete			Contract Manufacturer	Repackager / Relabeler
Company / Institution Nan	<u></u>		Establishment Registration Nu	mbor
			Latabianment Registration Nu	
Division Name (if applicat	le)	, ²⁰⁰	Phone Number (including area	a code)
Street Address			FAX Number (including area c	ode)
City			State / Province	ZIP Code Country
0,				
Contact Name		Contact Title		Contact E-mail Address
	Facility Establishment Identifier ((FEI) Number		
Original		, , ,	Manufacturer	Contract Sterilizer
Add Delete			Contract Manufacturer	Repackager / Relabeler
Company / Institution Nan	ne		Establishment Registration Nu	mber
Division Marra (Frank)			Dhana Numha- Cost atom	
Division Name (if applicat	ne)		Phone Number (including area	r coue)
Street Address			FAX Number (including area c	
				,
City	·····		State / Province	ZIP Code Country
	·	0	l	
Contact Name		Contact Title		Contact E-mail Address
FORM FDA 3514 (1/1	3)		Add	Continuation Page Page 4 of 5 Pages

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 307a396-8118f 55

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016

 Standards No.	Standards	Standards Title	Version	Date
11737-1	Organization	Sterilization of health care products - Microbiological methods - Part	2nd edition	
11/5/-1	ANSI/AAMI/ISO			04/01/2006
				04/01/2000
 Standards No.	Standards Organization	Standards Title	Version	Date
11737-2	ANSI/AAMI/ISO	Sterilization of health care products - Microbiological methods - Part	2nd edition	
		2		11/15/2009
 Standards No.	Standards	Standards Title	Version	Date
11137-2	Organization	Sterilization of health care products - Radiation - Part 2	2nd edition	
	ANSI/AAMI/ISO			03/15/2012
 Standards No.	Standards	Standards Title	Version	Date
11607-01	Organization	Packaging for terminally sterilized medical devices - Part 1	2006	
	ISO			01/01/2006
 Standards No.	Standards	Standards Title	Version	Date
see attached for additional standards	Organization	see attached for additional standards		
 Standards No.	Standards Organization	Standards Title	Version	Date
 Standards No.	Standards	Standards Title	Version	Date
	Organization			
 1	Please	include any additional standards to be cited on a separate pag	e.	J
		ion applies only to requirements of the Paperwork Reduction Act of 19 YOUR COMPLETED FORM TO THE PRA STAFF ADDRESS		
existing data source	es, gather and maintain	ormation is estimated to average 0.5 hour per response, including the n the data needed and complete and review the collection of informa- information collection, including suggestions for reducing this burden	tion. Send comm	
		Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff 1350 Piccard Drive, Room 400 Rockville, MD 20850		

FORM FDA 3514 (1/13)

Page 5 of 5 Pages

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 307a796-8218f 55

STANDARDS TO BE CITED:							
Item No.	Standard Number	Standards Organization	Standards Title	Version	Date		
6	10993-1	ISO	Biological evaluation of medical devices - Part 1: Evaluation and testing within a Risk Management Process	4 th edition	10/15/2009		
7	10993-5	ISO	Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity	3 rd edition	06/01/2009		
8	10993-10	ISO	Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization	3 rd edition	08/01/2010		
9	10993-11	ISO	Biological evaluation of medical devices – Part 11: Tests for systemic toxicity	2 nd edition	08/15/2006		
11	10555-1	IEN/ISO	Sterile, Single-use Intravascular Catheters – Part 1: General Requirements	1 st edition A2	05/15/2004		
12	1618	BS EN	Catheters other than intravascular catheters. Test methods for common properties	2001	2001		

3. 510(K) COVER LETTER

The cover letter is included on the following page.

(The remainder of this page is intentionally left blank.)



March 1, 2013

(b)(4)	
Attention: (b)(6)	
(b)(4)	

Re: Traditional 510(k) Notification - Broncus' LungPoint Tools

Dear((b)(6)

In accordance with 21 CFR 807 and FDA guidance entitled *Format for Traditional and Abbreviated 510(k)s* (August 2005), Broncus Medical, Inc, is submitting this Traditional 510(k) for the LungPoint Tools. The LungPoint Tools are Class II endoscopic tools used during bronchoscopic procedures.

The Traditional 510(k) is being submitted by Broncus to address the addition of two endoscopic accessories, referred to collectively as the LungPoint Tools. A previous submission was made to (b) for the LungPoint Tools as Special 510(k) to the LungPoint Software. Per our discussions, we are now submitting the LungPoint Tools as a separate Traditional 510(k).

The format of the submission is consistent with that required in the FDA guidance entitled *Format for Traditional and Abbreviated 510(k)s* (August 2005).

In accordance with 21 CFR 807.95, Broncus requests that **(b)** and FDA hold the contents of this submission as confidential. Broncus has taken precautions to protect the confidentiality of this information, and understands that the submission to the government of false information is prohibited by 18 USC 1001 and 221 USC 331(q).

Should	you	have	any	questions,	please	contact	me	at	(b)(4)	or	by	e-mail	at
(b)(4)													

Sincerely,

Consulting Regulatory Affairs Director

4. INDICATIONS FOR USE STATEMENT

The indication for use statement is included on the following page.

The intended use and indications of the LungPoint Tools, as described in its labeling, are the same as the intended use and indications for the predicate devices. **Section 12** provides a detailed comparison table.

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016 TRADITIONAL 510(k) NOTIFICATION

Indications for Use

510(k) Number (if known):	K
Device Name:	LungPoint Tools - LungPoint Dilation Balloon and LungPoint Sheath
Indications for Use:	The LungPoint Tools are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs.

Prescription Use X (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____ (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

Division Sign-Off Office of In Vitro Diagnostic Device Evaluation and Safety

510(k)_____

Page __ of ___

5. 510(K) SUMMARY

A 510(k) Summary of Safety and Effectiveness appears as **Appendix A** and is being submitted in accordance with the requirements of the Safe Medical Devices Act of 1990. This 510(k) Summary meets the requirements identified in 21 CFR 807.92.

(The remainder of this page is intentionally left blank.)

7. CLASS III SUMMARY AND CERTIFICATION

This section is not applicable since this 510(k) is for a Class II device and not for a device type classified into Class III.

(The remainder of this page is intentionally left blank.)

8. FINANCIAL CERTIFICATION OR DISCLOSURE STATEMENT

This section is not applicable since data from clinical studies is not used to support this submission.

(The remainder of this page is intentionally left blank.)

6. TRUTHFUL AND ACCURACY STATEMENT

TRUTHFUL AND ACCURATE STATEMENT

[As Required by 21 CFR 807.87(k)]

I certify that, in my capacity as Chief Operating Office of Broncus Medical, Inc., I believe to the best of my knowledge, that all data and information submitted in the premarket notification are truthful and accurate and that no material fact has been omitted.

Gary Kaplan Chief Operating Office

13

K______510(k) Number

9. DECLARATIONS OF CONFORMITY AND SUMMARY REPORTS

The development activities of the LungPoint Tools conformed to the following FDArecognized standard. FDA Forms 3654 for each standard listed are included in **Appendix I**.

Standard Number	Standards Title	Version					
Sterilization, Packaging and Shelf-life							
ANSI/AAMI/ISO 11737-1	Sterilization of health care products – Microbiological methods – Part 1: Determination of the population of microorganisms on product	2006					
ANSI/AAMI/ISO 11737-2	L = Dart 2: Tests of sterility performed in the validation of a						
ISO 11137-2	Sterilization of health care products – Radiation – Part 2:						
ISO 11607-1	Packaging for terminally sterilized medical devices – Part1: Requirements for materials, sterile barrier systems and packaging systems	2006					
Biocompatibility							
ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing within a Risk Management Process	2009					
ISO 10993-5	Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity	2009					
ISO 10993-10	Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization	2010					
ISO 10993-11	Biological evaluation of medical devices – Part 11: Tests for systemic toxicity	2009					
Bench Testing	· · · · ·	-					
IEN SO 10555-1	Sterile, Single-use Intravascular Catheters – Part 1: General Requirements	1995					
BS EN 1618	Catheters other than intravascular catheters - Test methods for common properties.	1997					

Table 9-1. List of Relevant Standards

11. DEVICE DESCRIPTION

The LungPoint Sheath and LungPoint Dilation Balloon are endoscopic tools, which are intended to be used as accessories to the LungPoint Software. They are described in more detail in Sections 11.1 and 11.2, respectively. The LungPoint Software may be used with these tools or other commercially available tools. However, these tools have been specifically designed to be used with the Software but can be used without the LungPoint Software, as well.

11.1. LungPoint Sheath

The LungPoint Sheath is an endoscopic tool that is designed to be used with bronchoscopes to provide a working channel through which endoscopic tools, such as needles, dilation balloons, or other endoscopic devices may be introduced to the targeted area within the respiratory organs.


The sheath $\binom{(b)(4)}{(b)(4)}$) is inserted through a standard bronchoscope. It is then advanced to a target under the guidance of the LungPoint Software. Placement of the sheath may require the use of other endoscopic tools, such as the dilation balloon described in the next section. Once placed, other endoscopic tools, such as needles, dilation balloons, or other endoscopic devices, may be introduced through the sheath to the targeted area within the respiratory organs. $\binom{(b)(4)}{(b)(4)}$

The LungPoint Sheath is provided as a sterile device for single use and is available in one size. A detailed drawing of the sheath is provided in **Appendix B**.

11.1.1. Materials

The materials used in the sheath are outlined below. All materials passed appropriate biocompatibility requirements as outlined in **Section 16**.

Table 11.1.1-1: LungPoint Sheath Materials List




Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016 TRADITIONAL 510(k) NOTIFICATION

Material Use	Material Description	Common Name
(4)		

Supporting Information

Material Safety Data Sheet (MSDS)

Certificate of Analysis



chemicals and their relative concentrations are included in the Material Safety Data Sheet (MSDS), as summarized in the table below.



Supporting information on the^{(b)(4)} is included in **Appendix C2** as outlined below.

	Supporting Information
	Material Safety Data Sheet (MSDS)
(D)(4	

11.1.2. Performance Specification

Performance specification of the device, acceptance criteria and rationale for the acceptance criteria for each specification to ensure safe and effective use are outlined below.

Performance Specification	Objective	Applicable Standards	Acceptance Criteria (confidence/reliability parameters)
Dimensional Testing	(b)(4)	NA – all criteria are per the product specifications	(b)(4)

Performance Specification	Objective	Applicable Standards	Acceptance Criteria (confidence/reliability parameters)
	(b)(4)		(b)(4)
Joints/ Tensile Test		Method: EN1618:1997 Acceptance Criteria: adapted from ISO 10555- 1:1995	
Simulated Use		NA – validated test fixture	

Performance Specification	Objective	Applicable Standards	Acceptance Criteria (confidence/reliability parameters)
Radiopacity	(b)(4)	NA	(b)(4)

11.2. LungPoint Dilation Balloon

The LungPoint Dilation Balloon is an endoscopic tool that is used to dilate the target lung

(0)(4)			
(b)(4)			

The LungPoint Dilation Balloon is provided as a sterile device for single use. When inflated (b)(4)

.

The materials used in the dilation balloon are outlined below.

Table 11.2-1: LungPoint Dilation Balloon Materials List

Material Use	Material Description	Trade Name
b)(4)		

Performance specification of the device and acceptance criteria for each specification to ensure safe and effective use are outlined below.

Test Title	Objective	Applicable Standards	Acceptance Criteria
Dimensional Testing	(b)(4)	NA – all criteria are per the product specifications	(b)(4)
Joints/ Tensile Test		Method: EN1618:1997	
		Acceptance Criteria: adapted from ISO 10555-1:1995	
Simulated Use		NA – validated test fixture	
Balloon burst pressure		NA	

Table 11.2-2: LungPoint Balloon Performance Specifications (D0210-007 B)

Broncus Medical, Inc.

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016 TRADITIONAL 510(k) NOTIFICATION

Test Title	Objective	Applicable Standards	Acceptance Criteria
Radiopacity	b)(4)	NA	(b)(4)
			·

10. EXECUTIVE SUMMARY


10.1. Background and Rationale

This Traditional 510(k) is being submitted by Broncus Medical, Inc to add two endoscopic tools – a sheath $\binom{(b)(4)}{4}$ and a dilation balloon - to the LungPoint Software System (version 3.0 of the software was previously cleared in combination with an aspiration needle called FlexNeedle under K112051).



10.2. Description of Device

The LungPoint Sheath and LungPoint Dilation Balloon are endoscopic tools used during



Additionally, both devices are standard bronchoscopic tools, which can be used without the LungPoint Software.

10.3. Device Comparison

The LungPoint Tools - LungPoint Sheath and Dilation Balloon - has the same intended use and technological characteristics as the predicate devices. Table 1 in Section 12 provides a detailed comparison of the LungPoint Sheath and Dilation Balloon and the predicate devices. Based on the information provided, Broncus believes that the LungPoint Tools do not introduce any new questions of safety and effectiveness and are substantially equivalent to the predicate devices.

10.4. Summary of Performance Testing

The design and safety of the LungPoint Sheath and LungPoint Dilation Balloon were verified by performing functional and performance testing. All tests were designed to subject the sheath and dilation balloon to stresses that exceed those which would be encountered during clinical use. Section 18 provides an overview of the performance testing conducted to support this premarket submission.

12. SUBSTANTIAL EQUIVALENCE DISCUSSION

12.1. Predicate Devices

The LungPoint Tools, which are designed to be used as accessories to the LungPoint Software System (with FlexNeedle) - cleared to market on October 12, 2011 (K112051), are commonly used endosopic tools. These tools can be used as stand-alone devices or with the use of software like the LungPoint Software or its predicate the inReach System.



The predicate devices identified are as follows:

510(k) Submitter	510(k) Number
Broncus Technologies, Inc*	K112051, cleared to market on October 12, 2011 (and K093423, K091160 and K090095, by reference)
SuperDimension	K110093, cleared to market on February 11, 2011 (and K071473, K092365 and K102604, by reference)
Boston Scientific Corporation	K023337, cleared to market on November 18, 2002
	Broncus Technologies, Inc* SuperDimension Boston Scientific

12.2. Comparison to the Predicate Devices

The following tables outline the similarities and differences between the LungPoint Tools and the predicate devices.

	able 12.2-1: Comparison Table - LungPoint Sneath				
		I-Logic inReach System	LungPoint Sheath		
		(K110093)			
(b)(4	Intended Use	No separate intended use exists for the guide sheath/catheter – it is sold as part of the inReach system and is used with the software and bronchoscopes to provide a pathway through which endoscopic tools, such as needles, dilation balloons, or other diagnostic/therapeutic devices may be introduced to the targeted area within the respiratory organs.	No separate intended use exists for the sheath – it is sold as part of the LungPoint Software System and is used to provide a pathway through which endoscopic tools, such as needles, dilation balloons, or other diagnostic/therapeutic devices may be introduced to the targeted area within the respiratory organs.		
	Radiopaque Markers on Sheath/Stylet	Yes	Yes		
(b)(4)	Sheath/Stylet				
	Delivered Through	Flexible bronchoscope with minimum working channel of by mm	Flexible bronchoscope with minimum working channel of .mm		
	Single Use	Yes	Yes		
	Sterilization Method)(4)			

Table 12.2-1: Comparison Table - LungPoint Sheath

	CRE Pulmonary Balloon Dilatation Catheter	LungPoint Dilation Balloon
Intended Use	(K023337) Used to endoscopically dilate strictures of the airway tree.	No separate intended use exists for the dilation balloon – as with the inReach System, it is sold as part of the Software System. The Dilation Balloon is used to endoscopically dilate the target tissue of the airway tree.
Balloon Size and Inflation (rated balloon pressure)*	(b) (4)	
Catheter Length		
Inflation Port		
Balloon Length*		
Balloon Inflation Method		
Delivery		
Single Use		
Sterilization Method		
b)(4)		

Table 12.2-2: Comparison Table - LungPoint Dilation Balloon

12.3. Substantial Equivalence Discussion

The LungPoint Sheath and Dilation Balloon, the subject devices of this premarket notification, are substantially equivalent to the predicate devices. As outlined above, the LungPoint Tools have the same intended use and technological characteristics as the predicate devices.

In making this determination, Broncus followed the FDA's Draft Guidance entitled The 510(k) Program: Evaluating Substantial Equivalence in Premarket Notifications [510(k)] (December 27, 2011), the detailed flowchart provided in the guidance (included at the end of this section) and 21 CFR 807.100.

Specifically the following questions were considered:

Decision 1: Is the predicate device legally marketed?

Yes all predicates are legally marketed. Specifically,

- the LungPoint Software System with FlexNeedle received clearance on October 12, 2011 (K112051),
- the inReach System received 510(k) clearance on February 11, 2011 (K110093), and
- the CRE Pulmonary Balloon received clearance on November 18, 2002 (K023337).

Decision 2: Do the devices have the same intended use?


The LungPoint Tools are accessories to the LungPoint Software System similar to (b)(4)

(b)(4)			

Decision 3: *Do the devices have the same technological characteristics?* Yes, the individual endoscopic tools provided with the software have the same technological characteristics as those currently on the market (refer to Table 12.2-1 and Table 12.2-2). Data supporting their performance, biocompatibility and sterility is included in the appropriate sections of this 510(k) submission.

The difference (b)(4)

Conclusions: Based on the analysis above, Broncus believes that the LungPoint Tools do not introduce any new questions of safety and effectiveness and are substantially equivalent to the predicate devices.



This Flowchart is not intended to be used as a 'stand-alone' document and should only be considered in conjunction with the accompanying text in this guidance.

Figure 1: 510(k) "Substantial Equivalence" Decision-Making Flowchart

13. PROPOSED LABELING

The labeling for the LungPoint Tools is provided in Appendix D. It includes the following:

- LungPoint Software and Tools Instructions for Use Addendum
- LungPoint Sheath Instructions for Use
- LungPoint Dilation Balloon Instructions for Use

All labeling was created in accordance with 21 CFR 801.

The 510(k) summaries for the predicate devices are included in Appendix E.

(The remainder of this page is intentionally left blank.)

15. BIOCOMPATIBILITY

Biocompatibility data is being submitted for the accessories being added to the LungPoint Software, specifically the LungPoint Sheath and LungPoint Dilation Balloon.

15.1. Overview

Biocompatibility testing was performed for both the LungPoint Sheath components



15.2. Materials in Contact with Patients

Table 15.2-1 and Table 15.2-2 provide an overview of the components of the LungPoint Sheath and LungPoint Dilation Balloon, respectively. The list identifies whether the devices

0)(4)		
		ł
	•	

Table 15.2-1: LungPoint Sheath Materials List

	Material Description	Material	Patient Contact (Yes/No)	
(b)	(4)			

	Material Description	Material	Patient Contact (Yes/No)
(b)	(4)		

Table 15.2-2: LungPoint Dilation Balloon Materials List

15.3. Biocompatibility Summary

The devices were evaluated for biocompatibility per ISO 10993-1:2009 requirements for Biological Evaluation of Medical Devices. Testing was performed (b) (b)(4) on post-sterilized devices per the contract testing facility's GLP protocols to meet the testing requirements summarized in Table 15.3-1.

Based on ISO 10993-1, the following tests were conducted.

	1 abit 13.3-1. Testing Requirements 150 10335-1		
	Test Article	ISO 10993-1, Biological evaluation of medical devices – Part 1: Evaluation and Testing	Biological Effect
(b)(4	4)		

Table 15.3-1: Testing Requirements ISO 10993-1

Table 15.3-2 summarizes the biocompatibility test results for the devices. Biocompatibility results show that both devices meet all test requirements. The complete test report is provided in **Appendix G.**

	Test Name	Test Description	Results	
(b)(4)				
	Test Report Refere	nce: (b)(4)		
	rest Report Refere			

Table 15.3-2: Testing Summary for LungPoint Tools

Based on the results of the biological testing, it can be concluded that the tissue contact with the devices does not result in adverse biological effects. Therefore, the LungPoint Tools are considered biocompatible, non-toxic, and safe for use as external communicating devices in limited contact with tissue.

14. STERILIZATION/SHELF LIFE

Sterilization data is being submitted for the LungPoint Sheath and LungPoint Dilation Balloon.

14.1. Overview

The LungPoint Sheath and LungPoint Dilation Balloon are both provided sterile for single

(b)(4)		
	Both products are ^{(b)(4)}) sterilized at (b)(4)	, a contract
(b)(4)	/`	(b)(4)

Upon receiving the devices back from the sterilization facility, the cartons containing the sterilized devices are removed from the corrugated box, and it is verified that incoming quantities match the number of originally sent devices. The certificate of irradiation is reviewed for appropriate dose delivery. The verification of sterilization acceptance is required before the devices are placed into Finished Goods.

The standards that apply to the LungPoint Tools sterilization processes include:

	Process	Standard	Standard Title
(b)	(4)		

Process	Standard	Standard Title	
(b)(4)			

The following sections describe the sterilization methods and validation for the Broncus devices.

14.2. Sterilization Process





14.6. Pyrogenicity

The LungPoint Tools (b)(4)

14.7. Shelf-Life

Broncus plans to introduce both the LungPoint Sheath and LungPoint Dilation Balloon with

)(4)

Testing of the LungPoint Tools includes the following performance specifications:

- Simulated Use Testing ((b)(4)
- Tensile (b)(4) testing)
- Balloon Fatigue and Burst Testing ((b)(4)

16. SOFTWARE

This section is not applicable as software is not included in this submission and the cleared LungPoint Software is not being modified to allow for the use of the LungPoint

(The remainder of this page is intentionally left blank.)

17. ELECTROMAGNETIC COMPATIBILITY/ELECTRICAL SAFETY

This section is not applicable for the device as there are no electrical components.

(The remainder of this page is intentionally left blank.)

CONFIDENTIAL

Broncus Medical, Inc.

19. PERFORMANCE TESTING – ANIMAL

Biocompatibility studies for the LungPoint Sheath and LungPoint Dilation Balloon were performed as outlined in Section 15. No additional animal studies were performed in support of this premarket submission given that no changes were made to the software as a result of the addition of the LungPoint Tools to the LungPoint Software.

(The remainder of this page is intentionally left blank.)

CONFIDENTIAL

Broncus Medical, Inc.

18. PERFORMANCE TESTING – BENCH

18.1. Overview

(b)(4)	

Requirements NOTE: Significance of each requirement has been modified form the FDA Guidance to remove references to stents	Applicable? Y/N	Rationale/Testing Conducted
)(4)		

Requirements		
NOTE: Significance of each requirement has been modified form the FDA	Applicable? Y/N	Rationale/Testing Conducted
(b)(4) Testing		

Requirements NOTE: Significance of each requirement has been modified form the FDA Guidance to remove references to stents	Applicable? Y/N	Rationale/Testing Conducted
(b)(4) Testing		

In addition to the testing outlined above, radiopacity of the markers on the LungPoint Tools was also confirmed using images of the catheter. Images of the markers were obtained during software developments testing and are included in **Appendix H5**.

Detailed test results of the LungPoint Dilation Balloon and LungPoint Sheath are outlined in the following sections. Appropriate standards, where applicable, were used to set acceptance criteria as outlined in **Section 11**.

18.2. LungPoint Dilation Balloon

The design and safety of the LungPoint Dilation Balloon was verified by performing functional and performance testing. Testing included:

		1.	Dimensional testing -(b)(b)(4) Testing (4)	
	(b)(4) Too	2.	Tensile (b)(4) (b _{testing} - (b)(4) Testing	
	(b)(4) Tes	SUN	9	
	:	3.	Burst pressure and fatigue $-(b)(4)$ Testing	
ting				

Table 18.2-1: Summary of Bench Testing for the Broncus Balloon Dilation Catheter

b)(4) Testing


CONFIDENTIAL Broncus Medical, Inc.

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 3(Page-4916f 55



18.3. LungPoint Sheath

The design and safety of the LungPoint Sheath was verified by performing functional and performance testing. Testing included:



b)(4) Testing

CONFIDENTIAL Broncus Medical, Inc.

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 3(Page-S116f 55

Table 18.3-1: Summary	of Bench Testing f	for the LungPoint S	heath
sting			

CONFIDENTIAL

Broncus Medical, Inc.

20. PERFORMANCE TESTING - CLINICAL

When compared to the predicate device, the LungPoint Software with the LungPoint Sheath and LungPoint Dilation Balloon does not incorporate any changes that impact safety or efficacy; therefore, as with the predicate device, a clinical study was not performed to support the premarket submission.

(The remainder of this page is intentionally left blank.)

Appendix A: 510(k) Summary

510(k) Summary

1. Date of Summary

March 1, 2013

2. 510(k) Applicant

Broncus Medical, Inc. 1400 N. Shoreline Blvd, Suite A8 Mountain View, California 94043 Phone: (650) 428-1600 FAX: (650) 428-1542

Contact Person:	Gary Kaplan
Phone:	(650) 428-1600
Fax:	(650) 428-1542
e-mail:	<u>gkaplan@broncus.com</u>

3. Device Overview

Trade Name:	LungPoint [™] Tools (LungPoint Sheath and LungPoint Dilation Balloon)
Common Name:	Sheath and Dilation Balloon
Classification Name:	Bronchoscope and Accessories 21 CFR 874.4680 Product Code EOQ

4. Predicate Device

The predicate devices identified are as follows:

Trade Name	510(k) Submitter	510(k) Number	
LungPoint Planning and Virtual Bronchoscopic Navigation (VBN) System (with FlexNeedle)	Broncus Technologies, Inc*	K112051, cleared to market on October 12, 2011 (and K093423, K091160 and K090095, by reference)	
inReach System	SuperDimension	K110093, cleared to market on February 11, 2011 (and K071473, K092365 and K102604, by reference)	
CRE Pulmonary Balloon Dilatation Catheter	Boston Scientific Corporation	K023337, cleared to market on November 18, 2002	
*NOTE: Broncus Technologies, Inc changed their name to Broncus Medical, Inc in June 2012.			

5. Device Description

The LungPoint Tools are endoscopic tools used during bronchoscopy procedures. The LungPoint Sheath is designed to be used with a bronchoscope to provide a working channel through which endoscopic tools, such as needles, dilation balloons, or other endoscopic devices may be introduced to the targeted lung tissue within the respiratory organs. It is advanced to a predefined target following guidance of the LungPoint Software. The LungPoint Dilation Balloon is used to dilate tissue of the bronchial tree and may be inserted through the sheath or directly through the working channel of the bronchoscope. Additionally, both devices could be used without the LungPoint Software.

6. Intended Use

The LungPoint Tools are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs.

7. Comparison to Predicate Device

The LungPoint Tools have the same intended use and technological characteristics as the predicate devices.

8. Performance Data

The design and safety of the LungPoint Tools were verified by performing functional and performance testing. All tests were designed to subject the sheath and dilation balloon to stresses that exceed those which would be encountered during clinical use.

9. Safety and Effectiveness

The LungPoint Tools labeling contains instructions for use and any necessary cautions and warnings, to provide for safe and effective use of the devices. The biocompatibility assessment of all patient contacting materials was performed in accordance with ISO 10993, *Biological Evaluation of Medical Devices*. In addition, the devices are sterilized using e-beam sterilization.

Appendix A

Appendix B: Device Drawings

Device drawing for the LungPoint Sheath (PN 12138) and LungPoint Balloon (PN 12139) are included.

Broncus Medical, Inc.

CONFIDENTIAL

Appendix C1: Supporting Information for ^{(b)(4)}

The following are included

- 1. MSDS (blue and green)
- 2. Certificate of Analysis (blue and green).

Appendix C2: Supporting Information for ^{(D)(4)}

The following are included

- 1. MSDS
- 2. Technical Data Sheet
- 3. Class VI Certification
- 4. Class Testing for Plastics

CONFIDENTIAL Broncus Medical, Inc.

Appendix D: Labeling

- 1. LungPoint Software and LungPoint Tools Instructions for Use Addendum (DRAFT)
 - 2. LungPoint Sheath Instructions for Use (DRAFT)
 - 3. LungPoint Dilation Balloon Instructions for Use (DRAFT)

CONFIDENTIAL Broncus Medical, Inc.

LungPoint Software with Endoscopic Tools Instructions for Use Addendum DRAFT

LungPoint Overview

The LungPoint Software is a device that guides a bronchoscope and endoscopic tools to a prespecified target in or adjacent to the bronchial tree by providing a path which is displayed on a 3D reconstruction of CT scan. The Software allows for visualization of the target in the lung tissue or in the bronchial tree; visualization of the interior of the bronchial tree; placement of catheters in the bronchial tree; and placement of markers into soft lung tissue to guide radiosurgery and thoracic surgery. The LungPoint Software is intended to be used with various commercially available endoscopic tools.

Consult the LungPoint Software User Manual for detailed information about LungPoint.

LungPoint Tools

Broncus Medical provides the following endoscopic tools (referred to as LungPoint Tools) for use during Virtual Bronchoscopic Navigation. A brief description of each tool follows the table. The LungPoint Tools are ebeam sterilized, supplied sterile, for single patient use.

LungPoint Tool	Part Number	Description
FlexNeedle	10005	18 gauge
LungPoint Dilation	TBD	Diameter x length: 4 x 6 mm
Balloon*		Catheter length: 155 cm
LungPoint Sheath*	TBD	ID/OD: 2.0/2.6 cm
		Catheter length: 100 cm

*NOTE: The LungPoint Sheath and Balloon are for use with the LungPoint Software only.

Caution: U.S. Federal law restricts the sale, distribution, and use of these devices to, or on the order of, a physician.

FlexNeedle

The FlexNeedle is an aspiration needle commonly used to collect specimens for histophatology or cytopathology. The needle has a coring needle tip in a protective sheath that is designed to fit down the 2-mm working channel of a standard flexible optical bronchoscope. The needle tip is welded to a coil tube that comprises a flexible shaft, allowing the needle to bend with the bronchoscope articulation and enabling access to hard-to-reach locations. When used together with the LungPoint Software, the needle can be safely guided to a prespecified targeted area within the respiratory organs in order to collect specimens.

LungPoint Dilation Balloon

The LungPoint Dilation Balloon is an endoscopic tool that is used to dilate the target lung tissue of the bronchial tree. The dilation balloon is inserted through the bronchoscope or sheath and is used to dilate lung tissue. The balloon may be used during interventional bronchoscopy procedures as determined by the physician performing the bronchoscopy.

LungPoint Sheath

The LungPoint Sheath is an endoscopic tool that is designed to be used with bronchoscopes to provide a working channel through which endoscopic tools, such as needles, dilation balloons, or other endoscopic devices may be introduced to the targeted area within the respiratory organs.

The sheath (with stylet, if necessary) is inserted through a standard bronchoscope. It is then advanced to a target under the guidance of the LungPoint Software. Placement of the sheath may require the use of other endoscopic tools, such as the dilation balloon described in the next section. Once placed, other endoscopic tools, such as needles, dilation balloons, or other endoscopic devices, may be introduced through the sheath to the targeted area within the respiratory organs. If a stylet is used, it must be removed prior to the introduction of the other tools.

WARNINGS AND CAUTIONS: LungPoint Tools

- The LungPoint Tools are supplied sterile and for single use only. Inspect the sealed sterile device package before opening. If the seal is broken, contents may not be sterile and could pose a risk of patient infection.
- If any damage or irregularity is found to the device after opening, DO NOT USE the device as this could result in harm to the patient.
- DO NOT USE the device past its "Use Before" (expiration) date.
- DO NOT USE the LungPoint Sheath with active devices.

Caution: U.S. Federal law restricts the sale, distribution, and use of these devices to, or on the order of, a physician.

Broncus Medical Inc.

LungPoint[™] Sheath

REF


Catalog number 10007

Instructions For Use (IFU)

Table of Contents

- **1** Device Description
- 2 Intended Use
- **3** Contraindications
- 4 Warnings / Precautions
- 5 **Possible Complications**
- 6 Preparation for Use
- 7 Operating Instructions
- 8 Storage

PN xxxxx_01 DRAFT



1 Device Description

The LungPoint Sheath is an endoscopic tool that is designed to be used with bronchoscopes to provide a working channel through which endoscopic tools, such as needles, dilation balloons, or other endoscopic devices may be introduced to the targeted area within the respiratory organs.

The sheath enables physicians to easily access targets and allows for multiple approaches to the preselected target. It allows for the repeated placement of endoscopic tools to a specified lesion(s) during one procedure. The sheath is provided with a stylet, which is used to minimize the amount of airway mucosa entering the sheath's lumen and/or to provide rigidity (i.e. pushability).

The sheath (with stylet, if necessary) is inserted through a standard bronchoscope. It is then advanced to a target under the guidance of the LungPoint Software. Placement of the sheath may require the use of other endoscopic tools. Once placed, other endoscopic tools, such as needles, dilation balloons, or other endoscopic devices, may be introduced through the sheath to the targeted area within the respiratory organs. If a stylet is used, it must be removed prior to the introduction of the other tools.

The LungPoint Sheath is comprised of a braid reinforced tubing, to resist kinking with articulation along with an incorporated stylet. The stylet when mated, provides for a rounded tip during insertion and pushability for navigation to the target. Removal of the stylet allows for standard 2.0mm working channel bronchoscopic accessories to be used though the lumen of the sheath. The tip of the Sheath is marked visually with black bands at the tip and in 1cm increments to 6cm to provide the user with an indication of the traveled depth during insertion. Also included are multiple radiopaque marker bands at the distal end @ the tip, 5mm, 10mm and 20mm to aid visualization of the sheath under fluoroscopy.

The LungPoint Sheath is designed to be used as a stand-alone device or with the LungPoint Software.

Broncus Medical Inc.

Product Configuration and Dimensions

The LungPoint Dilation Balloon (see illustration in Figure 1) has the following specifications:

Catalog Number	Working Length	Maximum Catheter Outer Diameter (OD)	Catheter Internal Diameter (ID)	Minimum accessory length for use through sheath
10007-1	900mm	2.65mm	2.0mm	965mm




Figure 1: LungPoint Sheath

2 Intended Use

The LungPoint Sheath is intended to provide a pathway through which endoscopic tools, such as needles, dilation balloons, or other diagnostic/therapeutic devices may be introduced to the targeted area within the respiratory organs.

3 Contraindications

There are no known contraindications.

R

Broncus Medical Inc.

(!) WARNINGS AND PRECAUTIONS

4 Warnings / Precautions

- The LungPoint Sheath should only be used by physicians thoroughly trained in bronchoscopy and in the particular technique and procedure to be performed, and familiar with the associated risks. These operating instructions must be read, understood, and followed.
- DO NOT USE the LungPoint Sheath with active devices.
- While advancing the sheath, DO NOT USE excessive force as this may advance the device in an uncontrolled manner that could result in harm to the patient's airway wall and cause bleeding.
- The LungPoint Sheath should only be used in a manner consistent with the Bronchoscope manufacturer's Instructions for Use. Use of the LungPoint Sheath in contradiction to the instructions may harm the device, bronchoscope or patient.

5 Possible Complications

- Bleeding
- Infection

6 Preparation for Use

6.1 Materials Required

- LungPoint Sheath
- Luer-lock aspiration syringe
- Flexible bronchoscope with a working channel of 2.8 mm or greater

6.2 Device Inspection and Preparation

6.2.1 Read all IFUs, package inserts, labels, and warnings for the device before beginning the procedure.

CAUTION: Inspect the sealed sterile device package before opening. If the seal is broken, contents may not be sterile and could pose a risk of patient infection.

6.2.2 Aseptically remove the LungPoint Sheath from the package and inspect for any damage, such as broken or crushed areas of the catheter shaft, sharp or protruding edges at the distal tip, or any kinks in the catheter shaft.

Broncus Medical Inc.

CAUTION: If any damage or irregularity is found, DO NOT USE the device as this could result in harm to the patient.

CAUTION: DO NOT USE the device past its "Use Before" (expiration) date.

6.2.3 Before inserting the device into the working channel of the bronchoscope, verify that the stylet is secured to the sheath via the luer lock connection.

7 Operating Instructions

7.1 Carefully advance the Sheath through the working channel of the bronchoscope using short (approximately 20 mm) increments, until the distal end of the catheter is within bronchoscopic view.

CAUTION: DO NOT force the catheter if resistance to insertion is encountered. Reduce the angulation of the bronchoscope until the device passes smoothly. This can result in (a) kinks to the device, (b) damage to the bronchoscope, or (c) harm to the patient such as punctures, hemorrhage or mucous membrane damage.

CAUTION: If any damage occurs, including a kink, DO NOT USE and discard.

7.2 Advance the catheter to the desired protrusion from the bronchoscope as required for navigation to the region of interest.

CAUTION: Do not angulate the bronchoscope abruptly while the device is extended from the distal end of the bronchoscope. This can result in kinks to the device or harm to the patient such as punctures, hemorrhage or mucous membrane damage.

7.3 For accessory use: remove the stylet by loosening the luer lock connection and withdraw the stylet. Insert the accessory device into the lumen.

CAUTION: If excessive resistance makes withdrawal difficult, reduce the angulation of the bronchoscope until the accessory can be withdrawn smoothly. Forcible withdrawal could damage the sheath, accessory or bronchoscope.

- 7.4 If using fluoroscopy visualize radiopaque accessory devices position relative to the tip marker band.
- 7.5 Withdraw the Sheath from the bronchoscope.
- 7.6 At the completion of the patient procedure, discard the device in the appropriate disposal container.

8 Storage

Store at controlled room temperature.



LEGAL NOTICE

© 2013 Broncus Medical, Inc. All rights reserved.

Broncus, the Broncus logo, and LungPoint are trademarks or registered trademarks of Broncus Medical, Inc.

Patents Pending.

REF	Catalog number
Ronly	Federal law (USA) restricts this device to sale by or on the order of a physician.
	Caution: consult instructions for use
STERILE R	Irradiation sterilized. Sterility guaranteed if package unopened and undamaged.
LOT	Serial number or batch code
	Use before
2	For single use
	Manufacturer

Graphic Symbol Legend for Medical Device Labeling



Manufacturer: Broncus Medical, Inc. 1400 N. Shoreline Blvd., Bldg. A, Suite 8 Mountain View, CA, 94043 USA (877) 428-1600

Broncus Medical Inc.

LungPoint[™] Dilation Balloon

REF

Catalog number 10008


Instructions For Use (IFU)

Table of Contents

- 1 Device Description
- 2 Intended Use
- 3 Contraindications
- 4 Warnings / Precautions
- 5 Possible Complications
- 6 Preparation for Use
- 7 Operating Instructions
- 8 Storage

PN xxxxx_01 DRAFT





1 Device Description

The LungPoint Dilation Balloon is an endoscopic tool that is used to dilate the target lung tissue of the bronchial tree. The dilation balloon is inserted through the bronchoscope or other endoscopic tool, such as a sheath, and is used to dilate lung tissue. The balloon may be used during interventional bronchoscopy procedures as determined by the physician performing the bronchoscopy.

The LungPoint Dilation Balloon is comprised of a catheter with a non-compliant 4mm OD x 6mm long inflatable balloon at the distal end and a y-connector with female luer fittings at the proximal end. The side port of the Y connector is used for inflation and deflation of the balloon with a standard indeflator capable of 12atm pressure (typically with a gage) and a male luer lock connector. The center port contains a nitinol wire which is fixed in place for pushability and kink resistance during tortuous bronchoscopic navigation. Included are two radiopaque marker bands and the distal and proximal ends of the balloon inflation length for visualization of the device position under fluoroscopy

The LungPoint Dilation Balloon is designed to be used as a stand-alone device or with the LungPoint Software.

Product Configuration and Dimensions

The LungPoint Dilation Balloon (see illustration in Figure 1) has the following specifications:

Catalog	Balloon Size	Rated Balloon	Maximum Catheter	Catheter
Number	(OD x length)	Pressure	Outer Diameter (OD)	Length/working length
10008-1	4mm x 6mm	12atm	1mm	1430mm/975mm

2 Intended Use

The LungPoint Dilation Balloon is intended to endoscopically dilate the target tissue of the airway tree.

LungPoint Dilation Balloon Instructions for Use Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796pendix D Page 11 of 15

Broncus Medical Inc.

3 Contraindications

There are no known contraindications.




Figure 1: Dilation Balloon (shown inflated)

(1) WARNINGS AND PRECAUTIONS

4 Warnings / Precautions

- The LungPoint Dilation Balloon should only be used by physicians thoroughly trained in bronchoscopy and in the particular technique and procedure to be performed, and familiar with the associated risks. These operating instructions must be read, understood, and followed.
- The LungPoint Dilation Balloon should only be used in a manner consistent with the Bronchoscope manufacturer's Instructions for Use. Use of the LungPoint Dilation Balloon in contradiction to the instructions may harm the device, bronchoscope or patient.

5 Possible Complications

- Bleeding
- Infection

6 Preparation for Use

- 5.1 Materials Required
 - LungPoint Dilation Balloon
 - Indeflator capable of 12atm
 - Flexible bronchoscope with a working channel of 2.0 mm or greater

LungPoint Dilation Balloon Instructions for Use Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796pendix D Page 12 of 15

K

Broncus Medical Inc.

- 5.2 Device Inspection and Preparation
 - 5.2.1 Read all IFUs, package inserts, labels, and warnings for the LungPoint Dilation Balloon before beginning the procedure.

CAUTION: Inspect the sealed sterile device package before opening. If the seal is broken, contents may not be sterile and could pose a risk of patient infection.

5.2.2 Aseptically remove the LungPoint Dilation Balloon from the package and inspect for any damage, such as broken or crushed areas of the catheter shaft, sharp or protruding edges at the distal tip, or any kinks in the catheter shaft.

CAUTION: If any damage or irregularity is found, **DO NOT USE** the device as this could result in harm to the patient.

CAUTION: DO NOT USE the device past its "Use Before" (expiration) date.

CAUTION: DO NOT place or manipulate the catheter in the airways unless under direct visualization with a bronchoscope as this may result in harm to the patient.

Broncus Medical Inc.

7 Operating Instructions

7.1 Carefully advance the balloon through the working channel of the bronchoscope using short (approximately 20 mm) increments, until the distal end of the catheter is within bronchoscopic view.

CAUTION: DO NOT force the catheter if resistance to insertion is encountered. Reduce the angulation of the bronchoscope until the device passes smoothly. This can result in (a) kinks to the device, (b) damage to the bronchoscope, or (c) harm to the patient such as punctures, hemorrhage or mucous membrane damage.

CAUTION: If any damage occurs, including a kink, DO NOT USE and discard.

7.2 Advance the catheter to the desired protrusion from the bronchoscope as required for navigation to the region of interest.

CAUTION: Do not angulate the bronchoscope abruptly while the device is extended from the distal end of the bronchoscope. This can result in kinks to the device or harm to the patient such as punctures, hemorrhage or mucous membrane damage.

- 7.3 Connect the indeflator via the luer lock connector and inflate to 10atm. Do not exceed the rated burst pressure of 12 atm.
- 7.4 If using fluoroscopy visualize balloon inflation location via the radiopaque marker bands and the distal and proximal ends of the balloon inflation length.
- 7.5 Deflate the balloon, and remove the device

CAUTION: If excessive resistance makes withdrawal difficult, reduce the angulation of the bronchoscope until the accessory can be withdrawn smoothly. Forcible withdrawal could damage the device or bronchoscope.

7.6 At the completion of the patient procedure, discard the device in the appropriate disposal container.

8 Storage

Store at controlled room temperature.

Broncus Medical Inc.

LEGAL NOTICE

© 2013 Broncus Medical, Inc. All rights reserved.

Broncus, the Broncus logo, and LungPoint are trademarks or registered trademarks of Broncus Medical, Inc.

Patents Pending.

REF	Catalog number
Ronly	Federal law (USA) restricts this device to sale by or on the order of a physician.
	Caution: consult instructions for use
STERILE R	Irradiation sterilized. Sterility guaranteed if package unopened and undamaged.
LOT	Serial number or batch code
	Use before
2	For single use
	Manufacturer

Graphic Symbol Legend for Medical Device Labeling



Broncus Medical Inc.

Manufacturer: Broncus Medical, Inc. 1400 N. Shoreline Blvd., Bldg. A, Suite 8 Mountain View, CA, 94043 USA (877) 428-1600

Appendix E: Predicate 510(k)Summaries

Summaries for the following devices are included,

- 1. LungPoint Software with FlexNeedle
- 2. inReach System
- 3. CRE Pulmonary Balloon Dilatation Catheter

CONFIDENTIAL Broncus Medical, Inc.

OCT. 13. 2011 2:48PM

NO. 9936 P. 1/3

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration 10903 New Hampshire Avenue Document Control Room – WO66-G609 Silver Spring, MD 20993-0002

Broncus Technologies, Inc. % Mr. Mark Job Responsible Third Party Official Regulatory Technology Services LLC 1394 25th Street NW BUFFALO MN 55313

OCT 1 2 2011

Re: K112051

Trade/Device Name: LungPoint[™] Planning and Virtual Brochoscopic Navigation (VBN) Software
Regulation Number: 21 CFR 892.2050
Regulation Name: Picture archiving and communications system
Regulatory Class: II
Product Code: LLZ
Dated: October 3, 2011
Received: October 6, 2011

Dear Mr. Job:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into class II (Special Controls), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016 Page 2

medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Parts 801 and 809), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbtanding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm for the CDRH's Office

of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <u>http://www.fda.gov/cdrh/industry/support/index.html</u>.

Sincerely Yours,

Mary Startel

Mary S. Pastel, Sc.D. Director Division of Radiological Devices Office of In Vitro Diagnostic Device Evaluation and Safety Center for Devices and Radiological Health

Enclosure

OCT. 13. 2011 2:49PM -

λ.

NO. 9936 P. 3/3

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016 *

. -

-

Indications for Use

510(k) Number (if known):	K112051
Device Name:	LungPoint™ Planning and Virtual Bronchoscopic Navigation (VBN) Software
Indications for Use:	Indicated for displaying images of the tracheobronchial tree to aid the physician in guiding endoscopic tools or catheters in the pulmonary tract and to enable marker placement within soft lung tissue. It does not make a diagnosis and is not an endoscopic tool. Not for pediatric use.

Prescription Use X (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____ (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

Division Sign-Off Office of In Vitro Diagnostic Device Evaluation and Safety

510(k) K112051

Page __ of ___

ŧ

Broncus Technologies, inc.

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-Appendix E Page 4 of 15

FEB 1 1 2011

510(k) SUMMARY superDimension® Marker Delivery Kit

superDimension, Ltd. 8 Hamenofim St. P.O. Box 2045 Herzliya 46120 Israel +972-(0)9-971-3700 +972-(0)9-971-3701 3004659744 Jonathan Kovach Vice President, Quality and Regulatory Affairs 763-210-4000 01/10/2011
P.O. Box 2045 Herzliya 46120 Israel +972-(0)9-971-3700 +972-(0)9-971-3701 3004659744 Jonathan Kovach Vice President, Quality and Regulatory Affairs 763-210-4000 01/10/2011
+972-(0)9-971-3701 3004659744 Jonathan Kovach Vice President, Quality and Regulatory Affairs 763-210-4000 01/10/2011 superDimension i Logic inReach System
3004659744 Jonathan Kovach Vice President, Quality and Regulatory Affairs 763-210-4000 01/10/2011 superDimension i Logic inReach System
Jonathan Kovach Vice President, Quality and Regulatory Affairs 763-210-4000 01/10/2011 superDimension i Logic inReach System
Vice President, Quality and Regulatory Affairs 763-210-4000 01/10/2011 superDimension i Logic inReach System
superDimension i Logic inReach System
superDimension inReach System superDimension/Bronchus
Bronchoscope
Computed tomography x-ray system
Radiology
21 CFR Part 892.1750
JAK
superDimension i Logic inReach System, K071473/K092365/K102604
Introduction of accessory Marker Delivery Kit for the placement of markers into soft lung tissue using the superDimension i Logic inReach System.
The superDimension i Logic inReach System is a device that guides a bronchoscope and endoscopic tools to a target in or adjacent to the bronchial tree on a path identified by CT scan. The superDimension i Logic inReach System also allows visualization of the target and the interior of the bronchial tree; placement of catheters in the bronchial tree; and placement of radiosurgical and dye markers into soft lung tissue to guide radiosurgery and thoracic surgery.
(

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-729peinetix E Page 5 of 15 •

dications for use	Indicated for displaying images of the tracheobronchial tree to aid the physician in guiding endoscopic tools or catheters in the pulmonary to to enable marker placement within soft lung tissue. It does not make diagnosis and is not an endoscopic tool. Not for pediatric use.				
Summary of the technological characteristics of the device compared to the predicate device					
Characteristic	Comparison to New Device				
Indications for Use	Same				
Anatomical Sites	Same				
Basic Principle	Same				
Localization Methodology	Same				
Visualization Principle and Method of Visualization of Path and Organs	Same				
Display Methods	Same				
Interventional Instrument	Same. Accessory Marker Delivery Kit is being added.				
Method of tracking location	Same				
Computer function and type	Same				
Energy Released into Body (Negligible)	Same				
Compatibility with other devices	Same. Accessory Marker Delivery Kit is being added.				
.ocompatibility	Same. Accessory Marker Delivery Kit is being added.				
Medical electrical equipment safety compliance	Same				
Registration Modes and Navigation Guidance	Same				
Software	Same				
Sterile Accessories	Same. Accessory Marker Delivery Kit is being added.				
Marker Placement	Same. Accessory Marker Delivery Kit is being added.				

SUMMARY OF NON-CLINICAL TESTS CONDUCTED FOR DETERMINATION OF SUBSTANTIAL EQUIVALENCE

Performance Test Summary-New Device

The superDimension Marker Delivery Kit and Instructions for Use were subjected to the superDimension design control process. Risk Management was performed to analyze the potential hazards associated with the changes. Appropriate design verification and validations were performed to assure the superDimension i Logic inReach System continues to be safe and effective for its intended use.

PERFORMANCE DATA

Comparative Performance Information Summary

Not Required to validate the changes to the superDimension System.

SUMMARY OF CLINICAL TESTS CONDUCTED FOR DETERMINATION OF SUBSTANTIAL EQUIVALENCE AND/OR OF CLINICAL INFORMATION

Clinical tests were not required to validate the changes to the superDimension System.

CONCLUSIONS DRAWN FROM NON-CLINICAL AND CLINICAL DATA

The superDimension i-Logic inReach System with the accessory Marker Delivery Kit is safe and effective for its intended use.

Page 3 of 3



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration 10903 New Hampshire Avenue Document Control Room – WO66-G609 Silver Spring, MD 20993-0002

superDimension, Ltd. % Mr Jonathan Kovach Vice President, Quality and Regulatory Affairs superDimension, Inc. 161 Cheshire Lane, Suite 100 MINNEAPOLIS MN 55441

FEB 1 1 2011

Re: K110093

Trade/Device Name: superDimension® Marker Delivery Kit Regulation Number: 21 CFR 892.1750 Regulation Name: Computed tomography x-ray system Regulatory Class: II Product Code: JAK Dated: January 10, 2011 Received: January 12, 2011

Dear Mr. Kovach:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into class II (Special Controls), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016

Page 2

medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Parts 801 and 809), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm</u> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address http://www.fda.gov/cdrh/industry/support/index.html.

Sincerely Yours,

Mary SVostel

Mary Pastel, ScD. Director Division of Radiological Devices Office of In Vitro Diagnostic Device Evaluation and Safety Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): |< 1/0093

Device Name: superDimension[®] Marker Delivery Kit

Indications for Use:

Indicated for displaying images of the tracheobronchial tree to aid the physician in guiding endoscopic tools or catheters in the pulmonary tract and to enable marker placement within soft lung tissue. It does not make a diagnosis and is not an endoscopic tool. Not for pediatric use.

Prescription Use _____ (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use ______(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

Page 1 of

(Division Sign-Off) Division of Radiological Devices Office of In Vitro Diagnostic Davice Evaluation and Safety

510K K110093

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-706perfetix E Page 10 of 15 Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016

NOV 1 8 2002

SECTION 14 510(K) SUMMARY

FOI RELEASABLE

Pursuant to §513(I)(3)(A) of the Food, Drug, and Cosmetic Act, Boston Scientific Corporation is required to submit with this Premarket Notification "...adequate summary of any information respecting safety and effectiveness or state that such information will be made available upon request of any person." Boston Scientific Corporation chooses to submit a summary of information respecting safety and effectiveness.

Date: Common/Usual Names: Trade/Proprietary Names:	October 4, 2002 Balloon Dilatation Catheter CRE [™] Pulmonary Balloon Dilatation Catheter		
Classification Name & Device Classification:	Class II		
	<u>Name</u> Bronchoscope & Acc.	<u>Number</u> 77 KTI	<u>21CFR Ref.</u> 874.4680
Device Panel/Branch:	Ear, Nose and Throat (ENT)		
Owner/Operator:	Boston Scientific Corporation One Boston Scientific Place Natick, MA 01760		
Contact Person:	Paige Sweeney Regulatory Affairs Special	st	

October 4, 2002 79 Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-749petridix E Page 11 of 15

Description of Devices

The CRETM Pulmonary Balloon Dilatation Catheter is used to access the airway tree via a bronchoscope for the purpose of dilating strictures. It consists of an inflatable balloon on a catheter shaft with lumens for inflation and guidewire.

Indications for Use

The CRE[™] Pulmonary Balloon Dilatation Catheter is intended to be used endoscopically to dilate strictures of the airway tree.

Descriptive and Technological Characteristics of Proposed and Predicate Devices

Boston Scientific Corporation believes that the CRE^{TM} Pulmonary Balloon Dilatation Catheter is substantially equivalent to the currently marketed CRE^{TM} Balloon Dilatation Catheter. The major components of these devices are the balloon, catheter shaft, and proximal hub. A thorough comparison of the descriptive characteristics between the proposed devices and the predicate devices show equivalence.

Conclusion

Boston Scientific Corporation has demonstrated that the CRE[™] Pulmonary Balloon Dilatation Catheter is substantially equivalent to the currently marketed CRE[™] Balloon Dilatation Catheter. Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016

DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

NOV 1 8 2002

Boston Scientific, Corp. c/o Paige Sweeney Regulatory Affairs Specialist Microvasive Endoscopy One Boston Scientific Place Natick, MA 01760

Re: K023337

Trade/Device Name: CRETM Pulmonary Balloon Dilation Catheter Regulation Number: 21 CFR 874.4680 Regulation Name: Bronchoscope (flexible or rigid) and accessories Regulatory Class: Class II Product Code: KTI Dated: October 4, 2002 Received: October 7, 2002

Dear Ms. Sweeney:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807): labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act): 21 CFR 1000-1050.

Page 2 – Paige Sweeney

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for <u>in vitro</u> diagnostic devices), please contact the Office of Compliance at (301) 594-4613. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address.

Sincerely yours,

A. Lalph forenthal

A. Ralph Rosenthal, M.D. Director Division of Ophthalmic and Ear, Nose and Throat Devices Office of Device Evaluation Center for Devices and Radiological Health

SECTION 3 INDICATIONS FOR USE

510(k) Number: To Be Determined $\angle 023337$

Device Name: - CRE[™] Pulmonary Balloon Dilatation Catheter

Indication for Use: The CRE[™] Pulmonary Balloon Dilatation Catheter is intended to be used endoscopically to dilate strictures of the airway tree.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE) Prescription Use _____ OR Over-The-Counter Use _____ (Per 21CFR 801.1091) (Optional Format 1-2-96)

(Division Sign-Off) Division of Ophthalmic Ear, Nose and Throat Devises

510(k) Number ______ K0 23337_____

510(k) Premarket Notification for CRE™ Pulmonary Balloon Dilatation Catheter Proprietary and Confidential Information of Boston Scientific Corporation

October 4, 2002

5

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-7Appetralix E Page 15 of 15 **Appendix F: Shelf-life Protocol**

CONFIDENTIAL Broncus Medical, Inc.

Appendix G: Biocompatibility Report QT-00804-R_A

CONFIDENTIAL Broncus Medical, Inc.

Appendix G Page 1 of 97 Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

Appendix H1: (b)(4) Testing

Broncus Medical, Inc.

CONFIDENTIAL

Appendix H2: (b)(4) Testing

Broncus Medical, Inc.

CONFIDENTIAL

Food and Drue STANDARDS DATA (To be filled i This report and the Summary Report Table are to be comp ences a national or international standard. A separate report TYPE OF 510(K) SUBMISSION	h and Human Services g Administration REPORT FOR 510(k)s n by applicant) pleted by the applicant when submitting a	510(k) th	at refer-			
STANDARD TITLE ¹ ISO 11737-1 Sterilization of medical devices - Microbiological m	ethods - Part 1: Determination of a population of	of microo	rganisms			
Please answer the following questions		Yes	No			
Is this standard recognized by FDA ² ?						
FDA Recognition number ³ # 14-326						
Was a third party laboratory responsible for testing conform in the 510(k)?	-	X				
Is a summary report ⁴ describing the extent of conformance 510(k)? If no, complete a summary report table.			X			
Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device?		X				
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).			X			
Does this standard include more than one option or selection of tests? If yes, report options selected in the summary report table.			X			
Were there any deviations or adaptations made in the use of the standard? If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) ⁵ ?						
Were deviations or adaptations made beyond what is specified in the FDA SIS?						
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			×			
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation Title of guidance:	of this 510k?					
¹ The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]	address of the test laboratory or certification body invo assessment to this standard. The summary report inc	ludes inforn				
² Authority [21 U.S.C. 360d], http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm	S.C. 360d], http://www.fda.gov/MedicalDevices/ nandGuidance/Standards/default.htm all standards utilized during the development of th 5 The supplemental information sheet (SIS) is addit is necessary before FDA recognizes the standard					
³ http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm ⁴ The summary report should include: any adaptations used to adapt to the	www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStanda	irds/search	.cfm			
⁴ The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and ⁶ The online search for CDRH Guidance Documents can be http://www.fda.gov/MedicalDevices/DeviceRegulationance GuidanceDocuments/default.htm						
FORM FDA 3654 (6/11) Questions? Contact FDA/CDRH/OCE/DID at CE	ne 1 ℃RH-FOISTATUS@fda.hhs.gov or 301-794		301) 443-6740 E			

Page 1 of 23

Recor	ds processed under FOIA Request # 2015 EXTENT OF STANDARD SUMMARY REPO	CONFORMANCE	02-01-201	16	• ••• ••••
	ion of medical devices - Microbiological metho	ods - Part 1			
	CONFORMANCE WITH STA	NDARD SECTIONS*			
SECTION NUMBER	SECTION TITLE		CONFORM	IANCE?	
1, 2, 3	Scope, Normative references, Terms and defin	nitions	🗶 Yes	🗌 No	N/A
TYPE OF DEVIATION OF NA	ROPTION SELECTED +		I		
DESCRIPTION					
JUSTIFICATION NA					
SECTION NUMBER	SECTION TITLE		CONFORM	ANCE?	
4	Quality managemetn system elements		🗶 Yes	No No	🗌 N/A
TYPE OF DEVIATION OF NA	ROPTION SELECTED *		<u>.</u>		
DESCRIPTION NA					
JUSTIFICATION NA					
SECTION NUMBER	SECTION TITLE		CONFORM	ANCE?	
5	Selection of product		🔀 Yes	🗌 No	N/A
TYPE OF DEVIATION OF NA	OPTION SELECTED *		L		
DESCRIPTION NA					
JUSTIFICATION NA					
explanation is needed described and adeque selected when follow report. More than on * Types of deviations of	all sections of the standard and indicate whether d under "justification." Some standards include of ately justified as appropriate for the subject devi ing a standard is required under "type of deviation e page may be necessary. can include an exclusion of a section in the stand S), a deviation to adapt the standard to the devi	options, so similar to deviations, th ice. Explanation of all deviations o on or option selected," "description dard, a deviation brought out by th	e option ch r descriptio " and "justi e FDA sup	osen nee on of optic ification" (eds to be ons on the
	Paperwork Reduction	Act Statement			
time for review completing and	g burden for this collection of information is est ing instructions, searching existing data sources reviewing the collection of information. Send co ollection of information, including suggestions	imated to average 1 hour per response of gathering and maintaining the date comments regarding this burden estimated	ata needed,	, and	
Department of Health and Human ServicesFood and Drug AdministrationOffice of Chief Information Officer1350 Piccard Drive, Room 400Rockville, MD 20850					

FORM FDA 3654 (6/11) Contact FDA/CDRH/OCE/DID at CBRH-FOISTATUS@fda.hhs.gov or 301-7 Apple hodix I Page 2 of 23 Form 3654 cont'd Standard Title: ISO 11737-1 Sterilization of medical devices - Microbiological methods - Part 1

(b)(4)

Appendix I: FDA Forms 3654

CONFIDENTIAL Broncus Medical, Inc.

Records processed under FOIA Request # 2015-7253, Regiessed by CDRH on 02-01-2015				
Department of Health Food and Drug STANDARDS DATA	n and Human Services g Administration REPORT FOR 510(k)s n by applicant)			
This report and the Summary Report Table are to be comp ences a national or international standard. A separate repor				
TYPE OF 510(K) SUBMISSION	Abbreviated		00.1974-1.197	
STANDARD TITLE ¹ ISO 11737-2 Sterilization of medical devices - Microbiological me	ethods - Part 2: Tests of sterility performed in the	ne definit	ion, valida	
Please answer the following questions		Yes	No	
Is this standard recognized by FDA ² ?		×		
FDA Recognition number ³		#_14-327		
Was a third party laboratory responsible for testing conform in the 510(k)?	•			
Is a summary report ⁴ describing the extent of conformance 510(k)? If no, complete a summary report table.			X	
Does the test data for this device demonstrate conformity to pertains to this device?				
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).			X	
Does this standard include more than one option or selection If yes, report options selected in the summary report table.	on of tests?	X		
Were there any deviations or adaptations made in the use of If yes, were deviations in accordance with the FDA supplem				
Were deviations or adaptations made beyond what is speci If yes, report these deviations or adaptations in the summar				
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			X .	
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation of Title of guidance:				
 ¹ The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication] ² Authority [21 U.S.C. 360d], http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm ³ http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm ⁴ The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and 	address of the test laboratory or certification body inv assessment to this standard. The summary report inc all standards utilized during the development of the d s The supplemental information sheet (SIS) is additional is necessary before FDA recognizes the standard. Fo www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStanda 6 The online search for CDRH Guidance Documents ca http://www.fda.gov/MedicalDevices/DeviceRegulation GuidanceDocuments/default.htm	ludes infor evice. al informatio ound at http ards/search an be found	mation on on which :// n.cfm d at	

FORM FDA 3654 (6/11) Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-740 EF Page 5 of 23

Records processed up	nder FOIA Reques	st # 2015-7253: Relea	ased by CDRH or	n 02-01-2016

	EXTENT OF STANDARD SUMMARY REPO					
STANDARD TITLE						
ISO 11737-2 Sterilization of medical devices - Microbiological methods - Part 2						
CONFORMANCE WITH STANDARD SECTIONS*						
SECTION NUMBER	SECTION TITLE		CONFORM	ANCE?		
1, 2, 3	Scope, Normative references, Terms and defin	nitions	🗙 Yes	🗌 No	□ N/A	
TYPE OF DEVIATION OF	R OPTION SELECTED +					
DESCRIPTION NA						
JUSTIFICATION NA						
SECTION NUMBER	SECTION TITLE		CONFORM	IANCE?	· · · · · · · · · · · · · · · · · · ·	
4	Quality management system elements		🗙 Yes	No	🗌 N/A	
TYPE OF DEVIATION OF	R OPTION SELECTED *					
DESCRIPTION						
JUSTIFICATION NA						
SECTION NUMBER	SECTION TITLE					
5	Selection of product		🗙 Yes	No No	□ N/A	
NA	R OPTION SELECTED *					
DESCRIPTION NA						
JUSTIFICATION NA						
explanation is neede described and adequ selected when follow report. More than or * Types of deviations of	t all sections of the standard and indicate wheth ad under "justification." Some standards include of uately justified as appropriate for the subject dev ving a standard is required under "type of deviation he page may be necessary. can include an exclusion of a section in the standard IS), a deviation to adapt the standard to the deviation	options, so similar to deviations, the ice. Explanation of all deviations o on or option selected," "description dard, a deviation brought out by th	e option ch r descriptio " and "justi	osen nee n of optic fication" d	eds to be ons on the	
		, ,	•			
time for review completing and	Paperwork Reduction g burden for this collection of information is est /ing instructions, searching existing data sources d reviewing the collection of information. Send ollection of information, including suggestions	timated to average 1 hour per response s, gathering and maintaining the date comments regarding this burden est	nta needed,	and		
Food Office 1350	tment of Health and Human Services and Drug Administration e of Chief Information Officer Piccard Drive, Room 400 ville, MD 20850	An agency may not conduct or spon required to respond to. a collection displays a currently valid OMB con	of informati	on unless		

FORM FDA 3654 (6(11)) Contact FDA/CDRH/OCE/DID at CORH-FOISTATUS@fda.hhs.gov or 301-7% pp Endix I Page 6 of 23 Form 3654 cont'd Standard Title: ISO 11737-2 Sterilization of medical devices - Microbiological methods - Part 2

(b)(4)

methods - Part 2

Form Approv	ed: OMR No	0910-0120-	Expiration	Date: 12/31/13
r onn Approv	EU. OMD NO.	. 0810-0120,	LApiration	Date. 12/01/10

Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)					
This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that refer- ences a national or international standard. A separate report is required for each standard referenced in the 510(k).					
TYPE OF 510(K) SUBMISSION		******			
Traditional Special	Abbreviated				
STANDARD TITLE ¹ 11137-2: Sterilization of health care products – Radiation – Part 2:	Establishing the sterilization dose				
Please answer the following questions		Yes	No		
Is this standard recognized by FDA ² ?		X			
FDA Recognition number ³		# <u>14-225</u>			
Was a third party laboratory responsible for testing conform in the 510(k)?		X			
Is a summary report ⁴ describing the extent of conformance 510(k)? If no, complete a summary report table.			\boxtimes		
Does the test data for this device demonstrate conformity to pertains to this device?	•	\boxtimes			
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).			X		
Does this standard include more than one option or selection of selection of the summary report table.	on of tests?		\boxtimes		
Were there any deviations or adaptations made in the use of If yes, were deviations in accordance with the FDA suppler					
Were deviations or adaptations made beyond what is speci If yes, report these deviations or adaptations in the summa					
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			X		
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation Title of guidance:					
¹ The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication] ² Authority [21 U.S.C. 360d], www.fda.gov/cdrh/stdsprog.html	certification body involved in conformance assessme standard. The summary report includes information of utilized during the development of the device.	on all stand			
 ³ http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/ search.cfm ⁴ The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or ⁵ The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/ search.cfm ⁶ The online search for CDRH Guidance Documents can be found at www.fda.gov/cdrh/guidance.html 					
FORM FDA 3654 (12/10) Pag	ge 1	PSC Graphic	s (301) 443-6740 EF		

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-7% poerformed in Page 8 of 23

Form Approved:	OMP No.	0010 0120-	Evniration	Data: 12/21/12
i unii Appiuveu.	UNID NO.	0910-0120,	LAPITATION	Uale. 12/01/10

Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)					
This report and the Summary Report Table are to be compl ences a national or international standard. A separate report					
TYPE OF 510(K) SUBMISSION	Abbreviated				
ISO 11607-1: Packaging for terminally sterilized medical devices -	Part 1: Requirements for materials, sterile barr	ier syste	ems		
Please answer the following questions		Yes	No		
Is this standard recognized by FDA ² ?		X			
FDA Recognition number ³		<u>14-355</u>			
Was a third party laboratory responsible for testing conformit in the 510(k)?		\boxtimes			
Is a summary report ⁴ describing the extent of conformance of 510(k)? If no, complete a summary report table.			X		
Does the test data for this device demonstrate conformity to pertains to this device?		\boxtimes			
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).			X		
Does this standard include more than one option or selection If yes, report options selected in the summary report table.	of tests?		\boxtimes		
Were there any deviations or adaptations made in the use of If yes, were deviations in accordance with the FDA suppleme					
Were deviations or adaptations made beyond what is specific If yes, report these deviations or adaptations in the summary					
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.	•		\boxtimes		
Is there an FDA guidance ⁶ that is associated with this standa If yes, was the guidance document followed in preparation of Title of guidance:					
 ¹ The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication] ² Authority [21 U.S.C. 360d], www.fda.gov/cdrh/stdsprog.html ³ http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/ search.cfm ⁴ The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or 	 certification body involved in conformance assessme standard. The summary report includes information or utilized during the development of the device. ⁵ The supplemental information sheet (SIS) is additionar which is necessary before FDA recognizes the standar http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfs search.cfm ⁶ The online search for CDRH Guidance Documents car www.fda.gov/cdrh/guidance.html 	n all stand al informa ard. Four Standards	tion id at /		

FORM FDA 3654 (12/10)

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-7%ppefridix I Page 9 of 23

	EXTENT OF STANDARI		
STANDARD TITLE ISO 11607-1: Packagii	ng for terminally sterilized medical devices - Pa	art 1: Requirements for materials,	sterile barrier systems
	CONFORMANCE WITH STA	ANDARD SECTIONS*	
SECTION NUMBER	SECTION TITLE		CONFORMANCE?
1, 2, 3, 4	Scope, Normative References, Terms and Det	initions, General requirements	X Yes No N/A
TYPE OF DEVIATION OF NA	ROPTION SELECTED *		
DESCRIPTION NA			
JUSTIFICATION NA			
SECTION NUMBER	SECTION TITLE		CONFORMANCE?
5	Materials and preformed sterile barrier system		Yes No N/A
TYPE OF DEVIATION OF NA	OPTION SELECTED *		L
DESCRIPTION NA			
JUSTIFICATION NA	· .		
SECTION NUMBER	SECTION TITLE	· · · · · · · · · · · · · · · · · · ·	CONFORMANCE?
6, 7	Design and development requirements for pac	kaging system, information	Yes No N/A
TYPE OF DEVIATION OF	ROPTION SELECTED *		
DESCRIPTION NA			
JUSTIFICATION NA			
explanation is neede described and adequ selected when follow report. More than on * Types of deviations of	all sections of the standard and indicate wheth d under "justification." Some standards include ately justified as appropriate for the subject dev ing a standard is required under "type of deviati e page may be necessary. can include an exclusion of a section in the stan S), a deviation to adapt the standard to the dev	options, so similar to deviations, th rice. Explanation of all deviations of on or option selected," "description dard, a deviation brought out by th	ne option chosen needs to be or description of options n" and "justification" on the ne FDA supplemental
	Paperwork Reduction	Act Statement	+
time for review completing and	g burden for this collection of information is es ing instructions, searching existing data source reviewing the collection of information. Send ollection of information, including suggestions	timated to average 1 hour per resp s, gathering and maintaining the d comments regarding this burden e	ata needed, and
Food a Office 1350 I	ment of Health and Human Services and Drug Administration of Chief Information Officer Piccard Drive, Room 400 ille, MD 20850	An agency may not conduct or spor required to respond to, a collection displays a currently valid OMB con	of information unless it
FORM FDA 3654 (12/1	D) Page 2	2	

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-7% performance for the Page 10 of 23

Form Approved: OMB No. 0910-0120; Expiration Date: 12/31/13

Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)					
This report and the Summary Report Table are to be comp ences a national or international standard. A separate report					
TYPE OF 510(K) SUBMISSION	44				
Traditional Special	Abbreviated				
STANDARD TITLE ¹ ISO 10993-1: Biological Evaluation of medical devices: Part 1: Eva	aluation and testing within a risk management	process			
Please answer the following questions		Yes	No		
Is this standard recognized by FDA ² ?		\mathbf{X}			
FDA Recognition number ³		# <u>2-179</u>			
Was a third party laboratory responsible for testing conformining the 510(k)?		\mathbf{X}			
Is a summary report ⁴ describing the extent of conformance 510(k)? If no, complete a summary report table.			\boxtimes		
Does the test data for this device demonstrate conformity to pertains to this device?	•	\times			
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).			X		
Does this standard include more than one option or selection If yes, report options selected in the summary report table.	n of tests?	\times			
Were there any deviations or adaptations made in the use or If yes, were deviations in accordance with the FDA supplem					
Were deviations or adaptations made beyond what is specifing the specific of the summary set of the summary			X		
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			X		
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation o Title of guidance:	f this 510k?				
 The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication] Authority [21 U.S.C. 360d], www.fda.gov/cdrh/stdsprog.html http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/ search.cfm The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods): choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or 	 certification body involved in conformance assessme standard. The summary report includes information of utilized during the development of the device. The supplemental information sheet (SIS) is addition which is necessary before FDA recognizes the stand http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cf search.cfm The online search for CDRH Guidance Documents c www.fda.gov/cdrh/guidance.html 	on all stand Ial informa Iard. Four Standards	tion nd at ;/		

FORM FDA 3654 (12/10)

Page 1

PSC Graphics (301) 443-6740 EF

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-7% ppferedix I Page 11 of 23

		DARD CONFORMANCE	
STANDARD TITLE ISO 10993-1: Biologi	cal Evaluation of medical devices: Part 1:	Evaluation and testing within a risk ma	anagement process
	CONFORMANCE WITH	STANDARD SECTIONS*	
SECTION NUMBER 1, 2, 3, 4	SECTION TITLE Scope, Normative References, Terms an	d Definitions, General Principles	CONFORMANCE?
TYPE OF DEVIATION O NA	R OPTION SELECTED *		
DESCRIPTION NA			
JUSTIFICATION NA			
SECTION NUMBER	SECTION TITLE Categorization of Medical Devices		
TYPE OF DEVIATION O	R OPTION SELECTED * are external communicating devices with I	imited (<74 hours) tissue contact por	$X Yes \qquad No \qquad N/A$
-			Section 5.2.2 and 5.5.
DESCRIPTION NA			
JUSTIFICATION NA			
SECTION NUMBER 6, 7	SECTION TITLE Biological Evaluation, Interpretation of B	iological Evaluation	CONFORMANCE?
	R OPTION SELECTED *		
NA			
DESCRIPTION NA			
JUSTIFICATION NA			
explanation is neede described and adeque selected when follow	at all sections of the standard and indicate we dounder "justification." Some standards inc uately justified as appropriate for the subject ving a standard is required under "type of d ne page may be necessary.	lude options, so similar to deviations, t t device. Explanation of all deviations	he option chosen needs to be or description of options
	can include an exclusion of a section in the IS), a deviation to adapt the standard to the		
	Paperwork Redu	iction Act Statement	• • • • • • • • • • • • • • • • • • •
time for review completing and	g burden for this collection of information ving instructions, searching existing data so I reviewing the collection of information. S ollection of information, including suggest	ources, gathering and maintaining the open comments regarding this burden	data needed. and
Food Office 1350	tment of Health and Human Services and Drug Administration e of Chief Information Officer Piccard Drive, Room 400 ville, MD 20850	An agency may not conduct or spo required to respond to, a collection displays a currently valid OMB co	n of information unless it
FORM FDA 3654 (12/1	0) P	age 2	······································

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-7% ppfeted ix I Page 12 of 23

r	A	OMD No.	0040 0400.	Europeanting.	Date: 10/01/10
rorm	Approvea:	ONB NO.	0910-0120;	Expiration	Date: 12/31/13

Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)					
This report and the Summary Report Table are to be comp ences a national or international standard. A separate report	· · · · · ·	• •	1		
TYPE OF 510(K) SUBMISSION					
Traditional Special	Abbreviated				
STANDARD TITLE ¹ ISO 10993-5: Biological Evaluation of medical devices: Part 5: Tes	st for in vitro cytotoxicity				
Please answer the following questions		Yes	No		
Is this standard recognized by FDA 2?		\times			
FDA Recognition number ³		# <u>2-153</u>			
Was a third party laboratory responsible for testing conformi in the 510(k)?	•	\boxtimes			
Is a summary report ⁴ describing the extent of conformance 510(k)? If no, complete a summary report table.			\boxtimes		
Does the test data for this device demonstrate conformity to pertains to this device?		\boxtimes			
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).		\boxtimes			
Does this standard include more than one option or selectio If yes, report options selected in the summary report table.	n of tests?	\boxtimes			
Were there any deviations or adaptations made in the use on If yes, were deviations in accordance with the FDA supplement					
Were deviations or adaptations made beyond what is specif If yes, report these deviations or adaptations in the summar					
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.	·····		\boxtimes		
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation o Title of guidance:		,			
 ¹ The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication] ² Authority [21 U.S.C. 360d], www.fda.gov/cdrh/stdsprog.html ³ http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/ search.cfm ⁴ The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or 	certification body involved in conformance assessm standard. The summary report includes information utilized during the development of the device. ⁵ The supplemental information sheet (SIS) is additio which is necessary before FDA recognizes the stan http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/c search.cfm ⁶ The online search for CDRH Guidance Documents www.fda.gov/cdrh/guidance.html	on all stan nal informa dard. Four cfStandards	ition nd at s/		

FORM FDA 3654 (12/10)

PSC Graphics (301) 443-6740 EF

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-7%ppehdix I Page 13 of 23

STANDARD TITLE	EXIENI OF SIA	NDARD CONFORMANCE Y REPORT TABLE			
	al Evaluation of medical devices: Part	5: Test for in vitro cytotoxicity			
	CONFORMANCE W	TH STANDARD SECTIONS*			
SECTION NUMBER 1, 2, 3	SECTION TITLE Scope, Normative References, Terms	and Definitions	CONFORMANCE?		
TYPE OF DEVIATION OF	ROPTION SELECTED *		L		
DESCRIPTION NA					
JUSTIFICATION NA					
SECTION NUMBER	SECTION TITLE		CONFORMANCE?		
4, 5, 6	Sample and control preparation, cell li	nes, culture medium	Yes No N/A		
TYPE OF DEVIATION OF	R OPTION SELECTED *				
DESCRIPTION					
JUSTIFICATION NA					
SECTION NUMBER	SECTION TITLE	and a duma Tast us ante Basult-	CONFORMANCE?		
	Preparation of cell stock culture, Test	procedures, l'est reports, Results	Yes No N/A		
TYPE OF DEVIATION OF Qualitative evaluation	was conducted (Section 8.5.1)				
DESCRIPTION NA					
JUSTIFICATION NA					
explanation is needed described and adequ selected when follow report. More than on	d under "justification." Some standards ately justified as appropriate for the sub ing a standard is required under "type o e page may be necessary.	te whether conformance is met. If a section include options, so similar to deviations, th oject device. Explanation of all deviations of f deviation or option selected," "description	ne option chosen needs to be or description of options n" and "justification" on the		
* Types of deviations can include an exclusion of a section in the standard, a deviation brought out by the FDA supplemental information sheet (SIS), a deviation to adapt the standard to the device, or any adaptation of a section.					
	Paperwork Re	eduction Act Statement	+		
time for review completing and	ing instructions. searching existing data	ion is estimated to average 1 hour per resp a sources, gathering and maintaining the d n. Send comments regarding this burden e gestions for reducing this burden to:	ata needed, and		
	tment of Health and Human Services				
Office 1350 F	and Drug Administration of Chief Information Officer Piccard Drive, Room 400 ille, MD 20850	An agency may not conduct or spor required to respond to, a collection displays a currently valid OMB con	of information unless it		
FORM FDA 3654 (12/1	0)	Page 2			

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-7% pperiodix I Page 14 of 23

E

Form Approved: OMB No. 0910-0120; Expiration Date: 12/31/13

Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)						
This report and the Summary Report Table are to be comple ences a national or international standard. A separate report						
TYPE OF 510(K) SUBMISSION						
🖂 Traditional 🔤 Special	Abbreviated					
STANDARD TITLE ¹ ISO 10993-10: Biological Evaluation of medical devices: Part 10: The second s	ests for irritation and skin sensitization					
Please answer the following questions		Yes	No			
Is this standard recognized by FDA ² ?		\boxtimes				
FDA Recognition number ³		<u></u> ¥2-174				
Was a third party laboratory responsible for testing conformit in the 510(k)?	•	\boxtimes				
Is a summary report ⁴ describing the extent of conformance of 510(k)? If no, complete a summary report table.		\boxtimes				
Does the test data for this device demonstrate conformity to t pertains to this device?		X				
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).		\boxtimes				
Does this standard include more than one option or selection If yes, report options selected in the summary report table.	of tests?	X				
Were there any deviations or adaptations made in the use of If yes, were deviations in accordance with the FDA suppleme	2					
Were deviations or adaptations made beyond what is specific If yes, report these deviations or adaptations in the summary						
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.	·····					
Is there an FDA guidance ⁶ that is associated with this standa If yes, was the guidance document followed in preparation of Title of guidance:						
 ¹ The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication] ² Authority [21 U.S.C. 360d], www.fda.gov/cdrh/stdsprog.html ³ http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/ search.cfm ⁴ The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or 	 certification body involved in conformance assessme standard. The summary report includes information of utilized during the development of the device. ⁵ The supplemental information sheet (SIS) is addition which is necessary before FDA recognizes the stand http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cf search.cfm ⁶ The online search for CDRH Guidance Documents of www.fda.gov/cdrh/guidance.html 	on all stand nal informa lard. Four Standards	tion nd at ;/			

FORM FDA 3654 (12/10)

PSC Graphics (301) 443-6740 EF

.

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-7% performance Page 15 of 23

		IDARD CONFORMANCE REPORT TABLE					
STANDARD TITLE ISO 10993-10: Biological Evaluation of medical devices: Part 10: Tests for irritation and skin sensitization							
CONFORMANCE WITH STANDARD SECTIONS*							
SECTION NUMBER	SECTION TITLE			CONFORMANCE?			
1, 2, 3, 4	Scope, Normative References, Terms a	nd Definitions, General Princi	ples	X Yes 🗌 No	N/A		
TYPE OF DEVIATION OF NA	ROPTION SELECTED *		1				
DESCRIPTION NA							
JUSTIFICATION NA							
SECTION NUMBER	SECTION TITLE			CONFORMANCE?			
5	Pretest consideration			🗙 Yes 🗌 No	🗌 N/A		
TYPE OF DEVIATION OF NA	OPTION SELECTED *		I				
DESCRIPTION NA							
JUSTIFICATION NA							
SECTION NUMBER	SECTION TITLE			CONFORMANCE?			
	Irritation Tests			Yes No	N/A		
TYPE OF DEVIATION OF Section 6.4 - Animal in	ROPTION SELECTED * tracutaneous (intradermal) reactivity test	was followed					
DESCRIPTION NA							
JUSTIFICATION NA							
explanation is needed described and adequa selected when followi report. More than one * Types of deviations c	all sections of the standard and indicate d under "justification." Some standards in ately justified as appropriate for the subje ng a standard is required under "type of o e page may be necessary. an include an exclusion of a section in the S), a deviation to adapt the standard to the	clude options, so similar to dev ect device. Explanation of all de deviation or option selected," "o e standard, a deviation brough	viations, the eviations or description' nt out by the	e option chosen ne description of opti and "justification"	eds to be ons on the		
1999-1997 - Tanana - T	Paperwork Red	uction Act Statement		¢			
time for reviewi completing and	burden for this collection of information ng instructions, searching existing data s reviewing the collection of information. llection of information, including sugges	ources, gathering and maintain Send comments regarding this	ning the da s burden es	ta needed, and	ŗ		
Food a Office 1350 P	ment of Health and Human Services nd Drug Administration of Chief Information Officer iccard Drive, Room 400 lle, MD 20850	An agency may not cond. required to respond to, a displays a currently valid	collection of	of information unless			
FORM FDA 3654 (12/10) F	Page 2					

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-7% ppetrolix I Page 16 of 23 Form 3654 cont'd Standard Title: ISO 10993-10: Biological Evaluation of medical devices - Part 10: Tests for irritation and skin sensitization



Records processed under FOIA Request # 2015-7253, ^F teileberereity/MBittein 89169189-59 Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)	ရာigtion Da	ate: 12/31/13
This report and the Summary Report Table are to be completed by the applicant when submitting a ences a national or international standard. A separate report is required for each standard referenced		
TYPE OF 510(K) SUBMISSION		
STANDARD TITLE ¹ ISO 10993-11: Biological evaluation of medical devices - Part 11: Tests for systemtic toxicity		
Please answer the following questions	Yes	No
Is this standard recognized by FDA ² ?	×	
FDA Recognition number ³	# <u>2-176</u>	
Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)?	×	
Is a summary report ⁴ describing the extent of conformance of the standard used included in the 510(k)? If no, complete a summary report table.		X
Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device?	X	
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).	X	
Does this standard include more than one option or selection of tests? If yes, report options selected in the summary report table.		X
Were there any deviations or adaptations made in the use of the standard? If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) ⁵ ?		
Were deviations or adaptations made beyond what is specified in the FDA SIS? If yes, report these deviations or adaptations in the summary report table.		
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.		
Is there an FDA guidance ⁶ that is associated with this standard? If yes, was the guidance document followed in preparation of this 510k? Title of guidance:		
 ¹ The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication] ² Authority [21 U.S.C. 360d], http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm ³ http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm ⁴ The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and 	cludes infor device. al informatio ound at http lards/search an be found	mation on on which ::// n.cfm d at

FORM FDA 3654 (6/11) Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-7% pendix I Page 18 of 23

Records	process	ed unde	r FOIA	Reque	est # 20	15-7253	: Relea	sed by	CDRH o	n 02-0 ⁻	1-2016
		with the March of the		10,0000-0010	· · · · · · · · · · · · · · · · · · ·			1.1.1.1.1.1.1.1.1.1	114.00 (14		

EXTEN	NT OF ST	FANDA	RD CON	IFORM	ANCE
	SUMMA	DV DEI	DODT T	ADI C	
				ADLE	

STANDARD TITLE ISO 10993-11: Biolog	gical evaluation of medical devices - Part 11: T	ests for systemtic toxicity			
SECTION NUMBER	SECTION TITLE		CONFORMANCE?		
1, 2, 3, 4	Scope, Normative references, Terms and def	initions, General considerations			
TYPE OF DEVIATION O	R OPTION SELECTED *				
DESCRIPTION NA					
JUSTIFICATION NA					
SECTION NUMBER	SECTION TITLE		CONFORMANCE?		
5	Acute Systemic Toxicity		Yes 🗌 No 🗌 N/A		
TYPE OF DEVIATION O	R OPTION SELECTED *		I		
DESCRIPTION NA					
JUSTIFICATION NA					
SECTION NUMBER	SECTION TITLE		CONFORMANCE?		
6	Repeated exposure systemic toxicity		Yes No N/A		
TYPE OF DEVIATION O	R OPTION SELECTED *		L		
DESCRIPTION NA					
JUSTIFICATION NA					
explanation is neede described and adequ selected when follow report. More than or * Types of deviations of	t all sections of the standard and indicate wheth d under "justification." Some standards include lately justified as appropriate for the subject de ring a standard is required under "type of deviat he page may be necessary. can include an exclusion of a section in the star IS), a deviation to adapt the standard to the dev	options, so similar to deviations, th vice. Explanation of all deviations of ion or option selected," "description ndard, a deviation brought out by th	e option chosen needs to be r description of options n" and "justification" on the e FDA supplemental		
Paperwork Reduction Act Statement					
time for review completing and	g burden for this collection of information is es ing instructions, searching existing data source reviewing the collection of information. Send ollection of information, including suggestions	es, gathering and maintaining the d comments regarding this burden e	ata needed, and		
Food a Office 1350 Rocky	tment of Health and Human Services and Drug Administration of Chief Information Officer Piccard Drive, Room 400 ville, MD 20850	An agency may not conduct or spon required to respond to. a collection displays a currently valid OMB con	of information unless it trol number.		
FORM FDA 3654 (6/11					

FORM FDA 365466(16) s? Contact FDA/CDRH/OCE/DID at CBRA-FOISTATUS@fda.hhs.gov or 301-796p81616dix I Page 19 of 23

Form Approved: OMB No. 0910-0120; Expiration Date: 12/31/13

Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)						
This report and the Summary Report Table are to be comp ences a national or international standard. A separate report						
TYPE OF 510(K) SUBMISSION						
Traditional Special	Abbreviated					
STANDARD TITLE ¹ ISO 10555-1: Sterile, Single-use Intravascular Catheters - Part 1: G	eneral Requirements (Amendment 2)					
Please answer the following questions		Yes	No			
Is this standard recognized by FDA 2?		\times				
FDA Recognition number ³		¥14-355				
Was a third party laboratory responsible for testing conformining the 510(k)?		\times				
Is a summary report ⁴ describing the extent of conformance of the standard used included in the 510(k)? If no, complete a summary report table.						
Does the test data for this device demonstrate conformity to pertains to this device?	-					
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).			X			
Does this standard include more than one option or selection If yes, report options selected in the summary report table.	n of tests?		X			
Were there any deviations or adaptations made in the use or If yes, were deviations in accordance with the FDA supplem						
Were deviations or adaptations made beyond what is specified of the summary of the summary set of the summar						
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			X			
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation o Title of guidance:						
 ¹ The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication] ² Authority [21 U.S.C. 360d], www.fda.gov/cdrh/stdsprog.html ³ http://www accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/ search.cfm ⁴ The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or 	 certification body involved in conformance assessme standard. The summary report includes information of utilized during the development of the device. ⁵ The supplemental information sheet (SIS) is addition which is necessary before FDA recognizes the stand http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cf search.cfm ⁶ The online search for CDRH Guidance Documents of www.fda.gov/cdrh/guidance.html 	on all stand al informa lard. Four Standards	tion id at i/			

FORM FDA 3654 (12/10)

Page 1

PSC Graphics (301) 443-6740 EF

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-7% performance Page 20 of 23

		NDARD CONFORMANCE	
STANDARD TITLE ISO 10555-1: Sterile, S	Single-use Intravascular Catheters - Par	rt 1: General Requirements (Amendment 2	2)
	CONFORMANCE W	/ITH STANDARD SECTIONS*	
SECTION NUMBER 1, 2, 3	SECTION TITLE Scope, Normative References, Defini	tions	CONFORMANCE?
TYPE OF DEVIATION OF	OPTION SELECTED *		I
DESCRIPTION NA			
JUSTIFICATION NA			
SECTION NUMBER	SECTION TITLE		CONFORMANCE?
	Requirements		X Yes No N/A
TYPE OF DEVIATION OF Sections 4.1, 4.2, and 4			
DESCRIPTION These sections refer to	general requirements, biocompatibility	and force at break.	to de d
JUSTIFICATION ISO 10555 is specific t sections were followed		ot required for bronchoscopic devices. The	erefore only the applicable
SECTION NUMBER	SECTION TITLE		CONFORMANCE?
	Description of nominal size, Informati	on to be supplied by the manufacturer	🗌 Yes 🗌 No 🔀 N/A
TYPE OF DEVIATION OF NA	ROPTION SELECTED *		
DESCRIPTION NA			
JUSTIFICATION ISO 10555 is specific to part of the standard is a		ot required for bronchoscopic devices. The	erefore compliance with only
 * For completeness list explanation is needed described and adequ selected when followi report. More than on * Types of deviations c 	all sections of the standard and indicat d under "justification." Some standards ately justified as appropriate for the sub ng a standard is required under "type o e page may be necessary. an include an exclusion of a section in	te whether conformance is met. If a sectior include options, so similar to deviations, th bject device. Explanation of all deviations o of deviation or option selected," "descriptior the standard, a deviation brought out by th the device, or any adaptation of a section.	e option chosen needs to be or description of options n" and "justification" on the ne FDA supplemental
	Danamuant D	aduation A at Statement	
time for review completing and	burden for this collection of informating instructions, searching existing data	eduction Act Statement ion is estimated to average 1 hour per resp a sources, gathering and maintaining the da n. Send comments regarding this burden es gestions for reducing this burden to:	ata needed, and
Food a Office 1350 P	ment of Health and Human Services nd Drug Administration of Chief Information Officer iccard Drive, Room 400 ille, MD 20850	An agency may not conduct or spon required to respond to, a collection displays a currently valid OMB con	of information unless it
FORM FDA 3654 (12/10))	Page 2	NAMA

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-7% pperfected ix I Page 21 of 23

Form Approved: OMB No. 0910-0120; Expiration Date: 12/31/13

Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)						
This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that refer- ences a national or international standard. A separate report is required for each standard referenced in the 510(k).						
TYPE OF 510(K) SUBMISSION						
🔀 Traditional 🔄 Special	Abbreviated					
STANDARD TITLE ¹ BS EN 1618: Catheters other than intravascular catheters - Test met	hods for common properties					
Please answer the following questions		Yes	No			
Is this standard recognized by FDA 2?			\boxtimes			
FDA Recognition number ³		<u><u></u>¥NA</u>				
Was a third party laboratory responsible for testing conformit in the 510(k)?			\mathbf{X}			
Is a summary report ⁴ describing the extent of conformance of the standard used included in the 510(k)? If no, complete a summary report table.						
Does the test data for this device demonstrate conformity to pertains to this device?	•	\boxtimes				
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).			X			
Does this standard include more than one option or selection If yes, report options selected in the summary report table.	of tests?		$\overline{\times}$			
Were there any deviations or adaptations made in the use of If yes, were deviations in accordance with the FDA suppleme						
Were deviations or adaptations made beyond what is specifi If yes, report these deviations or adaptations in the summary						
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.						
Is there an FDA guidance ⁶ that is associated with this stand If yes, was the guidance document followed in preparation of Title of guidance:	f this 510k?					
 ¹ The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication] ² Authority [21 U.S.C. 360d], www.fda.gov/cdrh/stdsprog.html ³ http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/ search.cfm ⁴ The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or 	certification body involved in conformance assessme standard. The summary report includes information of utilized during the development of the device. ⁵ The supplemental information sheet (SIS) is addition which is necessary before FDA recognizes the stand http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cf search.cfm ⁶ The online search for CDRH Guidance Documents of www.fda.gov/cdrh/guidance.html	on all stan al informa lard. Four Standards	tion nd at ;/			

FORM FDA 3654 (12/10)

Page 1

PSC Graphics (301) 443-6740 EF

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-7% performance Page 22 of 23

STANDARD TITLE		ANDARD CONFORMANCE RY REPORT TABLE	
BS EN 1618: Catheter	s other than intravascular catheters -	Test methods for common properties	
	CONFORMANCE	WITH STANDARD SECTIONS*	
SECTION NUMBER	SECTION TITLE Scope		CONFORMANCE?
TYPE OF DEVIATION OF	R OPTION SELECTED *		
DESCRIPTION NA			
JUSTIFICATION NA			
SECTION NUMBER 2, Annex A to F,ZA	SECTION TITLE Test Methods and Results		CONFORMANCE?
	R OPTION SELECTED * I to test methods for tensile property t	esting. Other appendices were not followed	ed
DESCRIPTION Annex B is the only an	nnex applicable to the LungPoint Too	ls.	
JUSTIFICATION The applicable annex v	was followed for tensile property test	ing.	
SECTION NUMBER	SECTION TITLE		
TYPE OF DEVIATION OF	R OPTION SELECTED *		Yes No N/A
DESCRIPTION			
JUSTIFICATION			
explanation is needed described and adequ selected when follow report. More than on	d under "justification." Some standard ately justified as appropriate for the s ing a standard is required under "type le page may be necessary.	ate whether conformance is met. If a secti is include options, so similar to deviations, ubject device. Explanation of all deviations of deviation or option selected," "description n the standard, a deviation brought out by	the option chosen needs to be s or description of options ion" and "justification" on the
		to the device, or any adaptation of a section	
	Paperwork	Reduction Act Statement	•
time for review completing and	ing instructions, searching existing date in the collection of information of the collection of the co	ation is estimated to average 1 hour per re ata sources, gathering and maintaining the ion. Send comments regarding this burden ggestions for reducing this burden to:	data needed, and
Food a Office 1350 F	tment of Health and Human Services and Drug Administration of Chief Information Officer Piccard Drive, Room 400 ille, MD 20850	An agency may not conduct or sp required to respond to, a collection displays a currently valid OMB c	on of information unless it
FORM FDA 3654 (12/10	0)	Page 2	

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-7%ppbehodix I Page 23 of 23

newborns hospitalized at this institution, including full access to the child, decision making authority for health care and other matters, and discharge of the child directly to them, and shall refer to them as parents of this child in any pertinent documents. Further, that the child's medical records shall be made available to Mahtab Fatemi or Bernard J.H. Verwer upon her or his request, without the need for further consent or authorization from any other entity;

- 7. The Ohio Division of Vital Records shall prepare the birth certificate consistent with this order, that is, the Division of Vital Records shall prepare the birth certificate of the child delivered by Melissa Wilcox on or about July 30, 2013 by referencing only the names and appropriate information of the judicially declared legal parents, Mahtab Fatemi and Bernard J.H. Verwer;
- Good cause exists to deem this case confidential and sealed from the public;
 a(b)(6)
 y, this Court orders that no person may examine the court file except the Parties and any attorney appoint by said Parties to act on their behalf and,
- 9. Such other and further relief as the nature of this cause may require.

Mahtab Fatemi,

Bernard J.H. Verwer

Melissa Wilcox	(b)(4)		
Michael Wilcox			




(b)(6)		
(b)(6)		



Broncus Medical, Inc. Broncus Medical, Inc. (650) 428-1600 1400 N. Shoreline Blvd. (650) 428-1542 fax Suite A8 Mountain View, CA 94043 www.broncus.com

(650) 428-1600

b)(4) Deficiencies



Broncus Medical, Inc. Page 2 of 15 CONFIDENTIAL Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

Attachment 1: CDRH Form 3564

Note: Submission of the inneed to submit device estimates	information entered in Section H does not affect the tablishment registration.	FDA Document Number (if kn	iown)
SECTION H	MANUFACTURING / PACKAGING /	STERILIZATION SITES RE	LATING TO A SUBMISSION
X Original	Facility Establishment Identifier (FEI) Number	X Manufacturer	Contract Sterilizer
Add Delete		Contract Manufacturer	Repackager / Relabeler
Company / Institution Nar	ne	Establishment Registration No	umber
Broncus Medical Inc		3007867778	
Division Name (if applical	ble)	Phone Number (including are	a code)
		650-428-1600	
Street Address		FAX Number (including area	code)
1400 N Shoreline Blvd,	A8	650-428-1542	

City State / Province ZIP Code Country Mountain View CA 94043 USA Contact Name Contact Title Contact E-mail Address John Magnasco VP, Operations and Quality Assurance jmagnasco@broncus.com Manufacturer Manufacturer Contract Sterilizer	ain View CA 94043 USA Name Contact Title Contact E-mail Address agnasco VP, Operations and Quality Assurance jmagnasco@broncus.com
Contact Name Contact Title Contact E-mail Address John Magnasco VP, Operations and Quality Assurance jmagnasco@broncus.com Imagnasco Facility Establishment Identifier (FEI) Number Manufacturer Imagnasco@broncus.com	Name Contact Title Contact E-mail Address agnasco VP, Operations and Quality Assurance jmagnasco@broncus.com inal Facility Establishment Identifier (FEI) Number Manufacturer X Contract Sterilizer
John Magnasco VP, Operations and Quality Assurance jmagnasco@broncus.com Imagnasco@broncus.com Imagnasco@broncus.com	inal Facility Establishment Identifier (FEI) Number
John Magnasco VP, Operations and Quality Assurance jmagnasco@broncus.com Imagnasco Facility Establishment Identifier (FEI) Number Imagnasco Imagnasco Manufacturer Imagnasco	inal Facility Establishment Identifier (FEI) Number
Imagination Imagination Imagination Imagination	inal Facility Establishment Identifier (FEI) Number
Contract Sterilizer	inal Manufacturer X Contract Sterilizer
Contract Sterilizer	inal Manufacturer X Contract Sterilizer

Original Add Delete	Facility Establishment Identifier (FEI) Number	Manufacturer Contract Manufacturer	Contract Sterilizer Repackager / Relabeler
Company / Institution Nan		Establishment Registration Nu	
Division Name (if annlinat		Dhong Number linduding area	a aadal

.

Division Name (if applicable)		Phone Number (including a	area code)	
Street Address		FAX Number (including are	a code)	
City		State / Province	ZIP Code	Country
Contact Name	Contact Title		Contact E-mail A	Address
FORM FDA 3514 (1/13)		Ac	Id Continuation Page	Page 4 of 5 Pages

Attachment 2: Indications for Use

Indications for Use

510(k) Number (if known):	K
Device Name:	LungPoint Tools - LungPoint Dilation Balloon and LungPoint Sheath
Indications for Use:	The LungPoint Tools are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs in a minimally invasive manner.

Prescription Use X (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____ (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page ____ of ____

Attachment 3: Instructions for Use LungPoint Sheath

Broncus Medical Inc.

LungPoint[™] Sheath

REF


Catalog number 10007

Instructions For Use (IFU)

Table of Contents

- **1** Device Description
- 2 Intended Use
- 3 Contraindications
- 4 Warnings / Precautions
- **5 Possible Complications**
- 6 Preparation for Use
- 7 Operating Instructions
- 8 Storage

PN xxxxx_01 DRAFT



1 Device Description

The LungPoint Sheath is an endoscopic tool that is designed to be used with bronchoscopes to provide a working channel through which endoscopic tools, such as needles, dilation balloons, or other endoscopic devices may be introduced to the targeted area within the respiratory organs.

The sheath enables physicians to easily access targets and allows for multiple approaches to the preselected target. It allows for the repeated placement of endoscopic tools to a specified lesion(s) during one procedure. The sheath is provided with a stylet, which is used to minimize the amount of airway mucosa entering the sheath's lumen and/or to provide rigidity (i.e. pushability).

The sheath (with stylet, if necessary) is inserted through a standard bronchoscope. It is then advanced to a target under the guidance of the LungPoint Software. Placement of the sheath may require the use of other endoscopic tools. Once placed, other endoscopic tools, such as needles, dilation balloons, or other endoscopic devices, may be introduced through the sheath to the targeted area within the respiratory organs. If a stylet is used, it must be removed prior to the introduction of the other tools.

The LungPoint Sheath is comprised of a braid reinforced tubing,to resist kinking with articulation along with an incorporated stylet. The stylet when mated, provides for a rounded tip during insertion and pushability for navigation to the target. Removal of the stylet allows for standard 2.0mm working channel bronchoscopic accessories to be used though the lumen of the sheath. The tip of the Sheath is marked visually with black bands at the tip and in <u>10mm</u> increments to <u>60mm</u> to provide the user with an indication of the traveled depth during insertion. Also included are multiple radiopaque marker bands at the distal end <u>at</u> the tip, 5mm, 10mm and 20mm to aid visualization of the sheath under fluoroscopy.

Broncus Medical Inc.

Product Configuration and Dimensions

The LungPoint Sheath (see illustration in Figure 1) has the following specifications:

Catalog Number	Working Length	Maximum Catheter Outer Diameter (OD)	Catheter Internal Diameter (ID)	Minimum accessory length for use through sheath
10007-1	900mm	2.65mm	2.0mm	965mm




Figure 1: LungPoint Sheath

2 Intended Use

The LungPoint Sheath is part of the LungPoint Tools, which are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs in a minimally invasive manner.

3 Contraindications

There are no known contraindications.

🕄 Broi

Broncus Medical Inc.

(!) WARNINGS AND PRECAUTIONS

4 Warnings / Precautions

- The LungPoint Sheath should only be used by physicians thoroughly trained in bronchoscopy and in the particular technique and procedure to be performed, and familiar with the associated risks. These operating instructions must be read, understood, and followed.
- <u>DO NOT activate active devices (e.g., electrosurgical devices) inside the lumen of</u> the sheath. The active portion of the device must be extended beyond the sheath prior to activation.
- While advancing the sheath, DO NOT USE excessive force as this may advance the device in an uncontrolled manner that could result in harm to the patient's airway wall and cause bleeding.
- The LungPoint Sheath should only be used in a manner consistent with the Bronchoscope manufacturer's Instructions for Use. Use of the LungPoint Sheath in contradiction to the instructions may harm the device, bronchoscope or patient.

5 Possible Complications

- Bleeding
- Infection

6 Preparation for Use

6.1 Materials Required

- LungPoint Sheath
- Luer-lock aspiration syringe
- Flexible bronchoscope with a working channel of 2.8 mm or greater

6.2 Device Inspection and Preparation

6.2.1 Read all IFUs, package inserts, labels, and warnings for the device before beginning the procedure.

CAUTION: Inspect the sealed sterile device package before opening. If the seal is broken, contents may not be sterile and could pose a risk of patient infection.

6.2.2 Aseptically remove the LungPoint Sheath from the package and inspect for any damage, such as broken or crushed areas of the catheter shaft, sharp or protruding edges at the distal tip, or any kinks in the catheter shaft.

Broncus Medical Inc.

CAUTION: If any damage or irregularity is found, DO NOT USE the device as this could result in harm to the patient.

CAUTION: DO NOT USE the device past its "Use Before" (expiration) date.

6.2.3 Before inserting the device into the working channel of the bronchoscope, verify that the stylet is secured to the sheath via the luer lock connection.

7 Operating Instructions

7.1 Carefully advance the Sheath through the working channel of the bronchoscope using short (approximately 20 mm) increments, until the distal end of the catheter is within bronchoscopic view.

CAUTION: DO NOT force the catheter if resistance to insertion is encountered. Reduce the angulation of the bronchoscope until the device passes smoothly. This can result in (a) kinks to the device, (b) damage to the bronchoscope, or (c) harm to the patient such as punctures, hemorrhage or mucous membrane damage.

CAUTION: If any damage occurs, including a kink, DO NOT USE and discard.

7.2 Advance the catheter to the desired protrusion from the bronchoscope as required for navigation to the region of interest.

CAUTION: Do not angulate the bronchoscope abruptly while the device is extended from the distal end of the bronchoscope. This can result in kinks to the device or harm to the patient such as punctures, hemorrhage or mucous membrane damage.

7.3 For accessory use: remove the stylet by loosening the luer lock connection and withdraw the stylet. Insert the accessory device into the lumen.

CAUTION: If excessive resistance makes withdrawal difficult, reduce the angulation of the bronchoscope until the accessory can be withdrawn smoothly. Forcible withdrawal could damage the sheath, accessory or bronchoscope.

- 7.4 If using fluoroscopy visualize radiopaque accessory devices position relative to the tip marker band.
- 7.5 Withdraw the Sheath from the bronchoscope.
- 7.6 At the completion of the patient procedure, <u>dispose of the device in</u> <u>accordance with applicable hospital, local, state and federal laws and</u> <u>regulations.</u>

8 Storage

Store at controlled room temperature.



LEGAL NOTICE

© 2013 Broncus Medical, Inc. All rights reserved.

Broncus, the Broncus logo, and LungPoint are trademarks or registered trademarks of Broncus Medical, Inc.

Patents Pending.

REF	Catalog number
Ronly	Federal law (USA) restricts this device to sale by or on the order of a physician.
	Caution: consult instructions for use
STERILE R	Irradiation sterilized. Sterility guaranteed if package unopened and undamaged.
LOT	Serial number or batch code
	Use before
2	For single use
	Manufacturer

Graphic Symbol Legend for Medical Device Labeling



Manufacturer: Broncus Medical, Inc. 1400 N. Shoreline Blvd., Bldg. A, Suite 8 Mountain View, CA, 94043 USA (877) 428-1600

Attachment 4: Instructions for Use LungPoint Dilation Balloon

Broncus Medical Inc.

LungPoint[™] Dilation Balloon

REF


Catalog number 10008

Instructions For Use (IFU)

Table of Contents

- 1 Device Description
- 2 Intended Use
- 3 Contraindications
- 4 Warnings / Precautions
- 5 Possible Complications
- 6 Preparation for Use
- 7 Operating Instructions
- 8 Storage





1 Device Description

The LungPoint Dilation Balloon is an endoscopic tool that is used to dilate the target lung tissue of the bronchial tree. The dilation balloon is inserted through the bronchoscope or other endoscopic tool, such as a sheath, and is used to dilate lung tissue. The balloon may be used during interventional bronchoscopy procedures as determined by the physician performing the bronchoscopy.

The LungPoint Dilation Balloon is comprised of a catheter with a non-compliant 4mm OD x 6mm long inflatable balloon at the distal end and a y-connector with female luer fittings at the proximal end. The side port of the Y connector is used for inflation and deflation of the balloon with a standard indeflator capable of 12atm pressure (typically with a gage) and a male luer lock connector. The center port contains a nitinol wire which is fixed in place for pushability and kink resistance during tortuous bronchoscopic navigation. Included are two radiopaque marker bands and the distal and proximal ends of the balloon inflation length for visualization of the device position under fluoroscopy

Product Configuration and Dimensions

The LungPoint Dilation Balloon (see illustration in Figure 1) has the following specifications:

Catalog	Balloon Size	Rated Balloon	Maximum Catheter	Catheter
Number	(OD x length)	Pressure	Outer Diameter (OD)	Length/working length
10008-1	4mm x 6mm	<u>20atm</u>	1mm	1430mm/975mm

3 Broncus Medical Inc.

2 Intended Use

The LungPoint Dilation Balloon is part of the LungPoint Tools, which are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs.

3 Contraindications

There are no known contraindications.




Figure 1: Dilation Balloon (shown inflated)

(1) WARNINGS AND PRECAUTIONS

4 Warnings / Precautions

- The LungPoint Dilation Balloon should only be used by physicians thoroughly trained in bronchoscopy and in the particular technique and procedure to be performed, and familiar with the associated risks. These operating instructions must be read, understood, and followed.
- The LungPoint Dilation Balloon should only be used in a manner consistent with the Bronchoscope manufacturer's Instructions for Use. Use of the LungPoint Dilation Balloon in contradiction to the instructions may harm the device, bronchoscope or patient.
- The rated burst pressure is 20 atmosphere, exceeding the rated burst pressure may result in harm to the patient.

B

Broncus Medical Inc.

5 Possible Complications

- Bleeding
- Infection

6 Preparation for Use

- 6.1 Materials Required
 - LungPoint Dilation Balloon
 - Indeflator capable of 12atm
 - · Flexible bronchoscope with a working channel of 2.0 mm or greater
- 6.2 Device Inspection and Preparation
 - 6.2.1 Read all IFUs, package inserts, labels, and warnings for the LungPoint Dilation Balloon before beginning the procedure.

CAUTION: Inspect the sealed sterile device package before opening. If the seal is broken, contents may not be sterile and could pose a risk of patient infection.

6.2.2 Aseptically remove the LungPoint Dilation Balloon from the package and inspect for any damage, such as broken or crushed areas of the catheter shaft, sharp or protruding edges at the distal tip, or any kinks in the catheter shaft.

CAUTION: If any damage or irregularity is found, **DO NOT USE** the device as this could result in harm to the patient.

CAUTION: DO NOT USE the device past its "Use Before" (expiration) date.

CAUTION: DO NOT place or manipulate the catheter in the airways unless under direct visualization with a bronchoscope <u>or using fluoroscopy</u> as this may result in harm to the patient.

Broncus Medical Inc.

7 Operating Instructions

7.1 Carefully advance the balloon through the working channel of the bronchoscope using short (approximately 20 mm) increments, until the distal end of the catheter is within bronchoscopic view.

CAUTION: DO NOT force the catheter if resistance to insertion is encountered. Reduce the angulation of the bronchoscope until the device passes smoothly. This can result in (a) kinks to the device, (b) damage to the bronchoscope, or (c) harm to the patient such as punctures, hemorrhage or mucous membrane damage.

CAUTION: If any damage occurs, including a kink, DO NOT USE and discard.

7.2 Advance the catheter to the desired protrusion from the bronchoscope as required for navigation to the region of interest.

CAUTION: Do not angulate the bronchoscope abruptly while the device is extended from the distal end of the bronchoscope. This can result in kinks to the device or harm to the patient such as punctures, hemorrhage or mucous membrane damage.

- 7.3 Connect the indeflator via the luer lock connector and inflate to 10atm. Do not exceed the rated burst pressure of 20 atm.
- 7.4 If using fluoroscopy visualize balloon inflation location via the radiopaque marker bands and the distal and proximal ends of the balloon inflation length.
- 7.5 Deflate the balloon, and remove the device

CAUTION: If excessive resistance makes withdrawal difficult, reduce the angulation of the bronchoscope until the accessory can be withdrawn smoothly. Forcible withdrawal could damage the device or bronchoscope.

7.6 At the completion of the patient procedure, <u>dispose of the device in</u> <u>accordance with applicable hospital, local, state and federal laws and</u> <u>regulations.</u>

8 Storage

Store at controlled room temperature.

Broncus Medical Inc.

LEGAL NOTICE

© 2013 Broncus Medical, Inc. All rights reserved.

Broncus, the Broncus logo, and LungPoint are trademarks or registered trademarks of Broncus Medical, Inc.

Patents Pending.

REF	Catalog number
Ronly	Federal law (USA) restricts this device to sale by or on the order of a physician.
	Caution: consult instructions for use
STERILE R	Irradiation sterilized. Sterility guaranteed if package unopened and undamaged.
LOT	Serial number or batch code
	Use before
2	For single use
	Manufacturer

Graphic Symbol Legend for Medical Device Labeling



Broncus Medical Inc.

Manufacturer: Broncus Medical, Inc. 1400 N. Shoreline Blvd., Bldg. A, Suite 8 Mountain View, CA, 94043 USA (877) 428-1600 Attachment 5: 510(k) Summary

510(k) Summary

1. Date of Summary

April 3, 2013

2. 510(k) Applicant

Broncus Medical, Inc. 1400 N. Shoreline Blvd, Suite A8 Mountain View, California 94043 Phone: (650) 428-1600 FAX: (650) 428-1542

Contact Person:	Gary Kaplan
Phone:	(650) 428-1600
Fax:	(650) 428-1542
e-mail:	<u>gkaplan@broncus.com</u>

3. Device Overview

Trade Name:	LungPoint [™] Tools (LungPoint Sheath and LungPoint Dilation Balloon)
Common Name:	Sheath and Dilation Balloon
Classification Name:	Bronchoscope and Accessories 21 CFR 874.4680 Product Code EOQ

4. Predicate Device

The predicate devices identified are as follows:

Trade Name	510(k) Submitter	510(k) Number	
LungPoint Planning and Virtual Bronchoscopic Navigation (VBN) System (with FlexNeedle)	Broncus Technologies, Inc*	K112051, cleared to market on October 12, 2011 (and K093423, K091160 and K090095, by reference)	
inReach System	SuperDimension	K110093, cleared to market on February 11, 2011 (and K071473, K092365 and K102604, by reference)	
CRE Pulmonary Balloon Dilatation Catheter	Boston Scientific Corporation	K023337, cleared to market on November 18, 2002	
*NOTE: Broncus Technologies, Inc changed their name to Broncus Medical, Inc in June 2012.			

5. Device Description

The LungPoint Tools are endoscopic tools used during bronchoscopy procedures. The LungPoint Sheath is designed to be used with a bronchoscope to provide a working channel through which endoscopic tools, such as needles, dilation balloons, or other endoscopic devices may be introduced to the targeted lung tissue within the respiratory organs. It is advanced to a predefined target following guidance of the LungPoint Software. The LungPoint Dilation Balloon is used to dilate tissue of the bronchial tree and may be inserted through the sheath or directly through the working channel of the bronchoscope. Additionally, both devices could be used without the LungPoint Software.

The materials used in the LungPoint Tools are commonly used medical grade materials and include platinum iridium markers.

6. Intended Use

The LungPoint Tools are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs in a minimally invasive manner.

7. Comparison to Predicate Device

The LungPoint Tools are commonly used endoscopic tools with the same technological characteristics as the predicate devices. The indications for use of the LungPoint Tools for use with the LungPoint Software all within the intended use of the predicate device, which is to enable treatment and diagnosis of lung cancer. The technological characteristics of the subject devices are the same as those of the predicate devices. Performance of the subject devices has been verified by use of accepted methods.

8. Performance Data

The design and safety of the LungPoint Tools were verified by performing functional and performance testing. All tests were designed to subject the sheath and dilation balloon to stresses that exceed those which would be encountered during clinical use. Testing included the following:

- Dimensional testing
- Joint/tensile test
- Simulated use
- Balloon fatigue/burst pressure
- Radiopacity.

9. Safety and Effectiveness

The LungPoint Tools labeling contains instructions for use and any necessary cautions and warnings, to provide for safe and effective use of the devices. The biocompatibility assessment of all patient contacting materials was performed in accordance with ISO 10993, *Biological Evaluation of Medical Devices*. In addition, the devices are sterilized using e-beam sterilization.

10. Conclusion

Based on the nonclinical testing the LungPoint Tools are as safe, as effective, and perform at least as safely and effectively as the predicate devices.

Appendix A

Attachment 6: Section 11 – Device Description

11. DEVICE DESCRIPTION

The LungPoint Sheath and LungPoint Dilation Balloon are endoscopic tools, which are intended to be used as accessories to the LungPoint Software. They are described in more detail in Sections 11.1 and 11.2, respectively. The LungPoint Software may be used with these tools or other commercially available tools. However, these tools have been specifically designed to be used with the Software but can be used without the LungPoint Software, as well.

11.1. LungPoint Sheath

The LungPoint Sheath is an endoscopic tool that is designed to be used with bronchoscopes to provide a working channel through which endoscopic tools, such as needles, dilation balloons, or other endoscopic devices may be introduced to the targeted area within the respiratory organs.

The sheath enables physicians to easily access targets and allows for multiple approaches to the preselected target. It allows for the repeated placement of endoscopic tools to a specified lesion(s) during one procedure. (b)(4)

<mark>b)(4)</mark>


The sheath (b)(4) is inserted through a standard bronchoscope. It is then advanced to a target under the guidance of the LungPoint Software. Placement of the sheath may require the use of other endoscopic tools, such as the dilation balloon described in the next section. Once placed, other endoscopic tools, such as needles, dilation balloons, or other endoscopic devices, may be introduced through the sheath to the targeted area within the respiratory organs. If a stylet is used, it must be removed prior to the introduction of the other tools.

The LungPoint Sheath is provided as a sterile device for single use and is available in one size. A detailed drawing of the sheath is provided in **Appendix B**.

11.1.1. Materials

The materials used in the sheath are outlined below. All materials passed appropriate biocompatibility requirements as outlined in **Section 16**.

Table 11.1.1-1: LungPoint Sheath Materials List



Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016 TRADITIONAL 510(k) NOTIFICATION

Material Use	Material Description	Common Name
(b)(4)		



Supporting information on the (b)(4) is included in Appendix C2 as outlined below.

11.1.2. Performance Specification

Performance specification of the device, acceptance criteria and rationale for the acceptance criteria for each specification to ensure safe and effective use are outlined below.

Performance Specification	Objective	Applicable Standards	Acceptance Criteria (confidence/reliability parameters)
Dimensional Testing	(b)(4)	NA – all criteria are per the product specifications	(b)(4)

 Table 11.1.2-1: LungPoint Sheath Performance Specifications (D0210-007 B)

Performance Specification	Objective	Applicable Standards	Acceptance Criteria (confidence/reliability parameters)	
	(b)(4)		(b)(4)	
Joints/ Tensile Test		Method: EN1618:1997		
		Acceptance Criteria: adapted from ISO 10555- 1:1995		
Simulated Use		NA – validated test fixture		

CONFIDENTIAL Broncus Medical, Inc. Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

Performance Specification	Objective	Applicable Standards	Acceptance Criteria (confidence/reliability parameters)
	(b)(4)		(b)(4)
Radiopacity		NA	

11.2. LungPoint Dilation Balloon

The LungPoint Dilation Balloon is an endoscopic tool that is used to dilate the target lung tissue of the bronchial tree. The dilation balloon is inserted through the bronchoscope or sheath and is used to dilate lung tissue. The balloon may be used (b)(4)

as determined by the physician performing the bronchoscopy.

The LungPoint Dilation Balloon is provided in one size, as outlined, and is a sterile device for single use.

	<u>Balloon Size</u> (OD x length)	<u>Balloon</u> <u>Pressure</u>	<u>Maximum Catheter</u> <u>Outer Diameter (OD)</u>	<u>Catheter</u> Length/working length	
(b)((4)				

When inflated at (b)(4) (diameter x length). A detailed drawing of the LungPoint Dilation Balloon is provided in **Appendix B**.

The materials used in the dilation balloon are outlined below.

Table 11.2-1: LungPoint Dilation Balloon Materials List

	Material Use	Material Description	Trade Name
(b)(4)			

CONFIDENTIAL Broncus Medical, Inc. Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

	Material Use	Material Description	Trade Name
(b)(4	1)		

Performance specification of the device and acceptance criteria for each specification to ensure safe and effective use are outlined below.


Test Title	Objective	Applicable Standards	Acceptance Criteria
Dimensional Testing	(b)(4)	NA – all criteria are per the product specifications	(b)(4)
Joints/ Tensile Test		Method: EN1618:1997 Acceptance Criteria: adapted from ISO 10555-1:1995	
Simulated Use		NA – validated test fixture	

 Table 11.2-2: LungPoint Balloon Performance Specifications (D0210-007 B)

CONFIDENTIAL Broncus Medical, Inc. Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

Test Title	Objective	Applicable Standards	Acceptance Criteria	
			(b)(4)	
Balloon burst pressure	(b)(4)	NA		
Radiopacity	(b)(4) (b)(4)	NA		

Attachment 7: Section 12 – Substantial Equivalence Discussion



This Flowchart is not intended to be used as a 'stand-alone' document and should only be considered in conjunction with the accompanying text in this guidance.

Figure 1: 510(k) "Substantial Equivalence" Decision-Making Flowchart

CONFIDENTIAL Broncus Medical, Inc. **Attachment 8: Pouch and Package Labels**



Broncus Medical Inc.

CONFIDENTIAL

All in ormation contained in this document is con idential and is the sole property o Broncus Medical Any reproduction in part or whole without the written permission o Broncus Medical is prohibited



Sheath

NOTE THAT THIS DOCUMENT IS ONLY FOR REVISION HISTORY OF 12198 IMAGE FILES

REVISION HISTORY – 12198 (CONTENT)

Rev	DCN #	Change Description	Release Date
		Initial release	11/14/12

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016



Broncus Medical Inc.



All in ormation contained in this document is con idential and is the sole property o Broncus Medical Any reproduction in part or whole without the written permission o Broncus Medical is prohibited



Label, Pouch, LungPoint Sheath



Broncus Medical, Inc.

1400 N. Shoreline Blvd., Bldg. A, Suite 8 Mountain View, CA 94043 USA 877-428-1600 www.broncus.com

LungPoint[™] Sheath



Sheath OD: 2.6 mm, Sheath ID: 2.0 mm Working Length: 900 mm



Ronly

For single use. Sterility guaranteed STERILE if package unopened and undamaged.



Rev.A

12198

P/N:



Broncus Medical, Inc. 1400 N. Shoreline Blvd., Bldg. A, Suite 8

Mountain View, CA 94043 USA 877-428-1600 © 2012 Broncus Medical, Inc. All rights reserved. Patent Pending. Broncus, B (logo), and LungPoint are trademarks or registered trademarks of Broncus Medical, Inc.



See package insert for explanation of symbols.


REF	10007		
LOT			
\blacksquare			
REF	10007		
LOT			
\square			

	Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016		
ß	Broncus Medical Inc.	Part No: #: <mark>(b)(4)</mark> Rev: <mark>(</mark> b	RELEASED
CONFIL	DENTIAL		
All in ormation contained in this document is con idential and is the sole property o Broncus Medical Any reproduction in part or whole without the written permission o Broncus Medical is prohibited		Label, Poue LungPoint Balloon	•

NOTE THAT THIS DOCUMENT IS ONLY FOR REVISION HISTORY OF 12200 IMAGE FILES.

REVISION HISTORY – 12200 (CONTENT)

Rev	DCN #	Change Description	Release Date
((b)	Initial release	11/14/12



Attachment 9: LungPoint Tools Addendum

LungPoint Software with Endoscopic Tools Instructions for Use Addendum DRAFT

LungPoint Overview

The LungPoint Software is a device that guides a bronchoscope and endoscopic tools to a prespecified target in or adjacent to the bronchial tree by providing a path which is displayed on a 3D reconstruction of CT scan. The Software allows for visualization of the target in the lung tissue or in the bronchial tree; visualization of the interior of the bronchial tree; placement of catheters in the bronchial tree; and placement of markers into soft lung tissue to guide radiosurgery and thoracic surgery. The LungPoint Software is intended to be used with various commercially available endoscopic tools.

Consult the LungPoint Software User Manual for detailed information about LungPoint.

LungPoint Tools

Broncus Medical provides the following endoscopic tools (referred to as LungPoint Tools) for use during Virtual Bronchoscopic Navigation. A brief description of each tool follows the table. The LungPoint Tools are ebeam sterilized, supplied sterile, for single patient use.

LungPoint Tool	Part Number	Description
FlexNeedle	10005	18 gauge
LungPoint Dilation	TBD	Diameter x length: 4 x 6 mm
Balloon*		Catheter length: <u>143 cm</u>
LungPoint Sheath*	TBD	ID/OD: 2.0/2.6 cm
		Catheter length: <u>975 mm</u>

*NOTE: The LungPoint Sheath and Balloon are for use with the LungPoint Software only.

FlexNeedle

The FlexNeedle is an aspiration needle commonly used to collect specimens for histophatology or cytopathology. The needle has a coring needle tip in a protective sheath that is designed to fit down the 2-mm working channel of a standard flexible optical bronchoscope. The needle tip is welded to a coil tube that comprises a flexible shaft, allowing the needle to bend with the bronchoscope articulation and enabling access to hard-to-reach locations. When used together with the LungPoint Software, the needle can be safely guided to a prespecified targeted area within the respiratory organs in order to collect specimens.

LungPoint Dilation Balloon

The LungPoint Dilation Balloon is an endoscopic tool that is used to dilate the target lung tissue of the bronchial tree. The dilation balloon is inserted through the bronchoscope or sheath and is used to dilate lung tissue. The balloon may be used during interventional bronchoscopy procedures as determined by the physician performing the bronchoscopy.

LungPoint Sheath

The LungPoint Sheath is an endoscopic tool that is designed to be used with bronchoscopes to provide a working channel through which endoscopic tools, such as needles, dilation balloons, or other endoscopic devices may be introduced to the targeted area within the respiratory organs.

The sheath (with stylet, if necessary) is inserted through a standard bronchoscope. It is then advanced to a target under the guidance of the LungPoint Software. Placement of the sheath may require the use of other endoscopic tools, such as the dilation balloon described in the <u>previous</u> section. Once placed, other endoscopic tools, such as needles, dilation balloons, or other endoscopic devices, may be introduced through the sheath to the targeted area within the respiratory organs. If a stylet is used, it must be removed prior to the introduction of the other tools.

WARNINGS AND CAUTIONS: LungPoint Tools

- The LungPoint Tools are supplied sterile and for single use only. Inspect the sealed sterile device package before opening. If the seal is broken, contents may not be sterile and could pose a risk of patient infection.
- If any damage or irregularity is found to the device after opening, DO NOT USE the device as this could result in harm to the patient.
- DO NOT USE the device past its "Use Before" (expiration) date.
- <u>Consult the individual instructions for use of the LungPoint Tools for additional</u> <u>warnings and cautions.</u>

Caution: Federal Law restricts this device to sale by or on the order of a physician DRAFT 25 March 2013 Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118 **Attachment 10: Updated Shelf-life Protocol**

Attachment 11: Section 18 – Bench Testing



b)(4) Schematic Drawings

b)(4) Schematic Drawings

Attachment 13: Simulated Use Testing Flowchart (b)(4) Testing



Flowchart 1

Attachment 14: Predicate Device Labeling

Included in this attachment are the following,

- 1. Marketing Literature for the inReach System Components and Tools
- K110093 inReach System
 K102604 inReach System
- 4. Product Label for CRE Balloon

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016 Components & Tools

superDimension leads the Electromagnetic Navigation Bronchoscopy[™] (ENB) market with its i·Logic System. Designed to enable physicians to reach distant lung lesions in a minimally invasive manner, the i·Logic System consists of intuitive software and proprietary hardware.

i-Logic Components


Planning & Navigation Software

- A disposable Extended Working Channel (EWC) that extends beyond the reach of the bronchoscope and becomes a channel for endobronchial tools to distal locations in the lungs and mediastinum
- A disposable Locatable Guide (LG) that contains a **location sensor** at its distal tip and allows 360 degree steerability through the bronchial tree
- Planning and navigation software, providing physicians with a reconstruction of the bronchial airways of the lungs as well as tracking and localization through Electromagnetic Navigation Bronchoscopy[™] (ENB).



Tools

The superDimension superTrax[™] tools are the only tools specifically designed for use with the i·Logic System. These sterile single use tools feature flexible shaft designs for maximum trackability through distal airways.



Ordering Information for Hogic[™] Tools Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016

DESCRIPTION	DIAMETER (mm)	WORKING LENGTH (cm)	REQUIRED WORKING CHANNEL (mm)	UNIT
Biopsy Forceps	1.7	110	2.0	5 / box
Cytology Brushes	1.7	120	2.0	10 / box
21g Aspirating Needles	1.8	130	2.0	5 / box

Ordering Information for i-Logic[™] System and Components

DESCRIPTION
i-Logic™ System
i-Logic™ All-Inclusive Olympus* Procedure Kit
i·Logic™ All-Inclusive Pentax® Procedure Kit
i-Logic [™] All-Inclusive Fujinon [®] Procedure Kit
i-Logic [™] Steerable Navigation Catheter (LG)
i-Logic™ Guide Catheter (EWC)
Fujinon® Bronchoscope Adaptor
Olympus®Bronchoscope Adaptor
Pentax [®] Bronchoscope Adaptor
Bronchoscope Clip
Patient Sensor Patches (box of 60)

*i·Logic is a trademark of superDimension, Ltd. | Disclaimer

FEB 1 7 2011

.

510(k) SUMMARY superDimension® Marker Delivery Kit

Name	superDimension, Ltd.		
Address	8 Hamenofim St. P.O. Box 2045 Herzliya 46120 Israel		
Phone number	+972-(0)9-971-3700 +972-(0)9-971-3701		
Fax number			
Establishment Registration Number	3004659744		
Name of contact person	Jonathan Kovach Vice President, Quality and Regulatory Affairs 763-210-4000		
Date prepared	01/10/2011		
Name of device			
Trade or proprietary name	superDimension i Logic inReach System superDimension inReach System superDimension/Bronchus		
Common or usual name	Bronchoscope		
Classification name	Computed tomography x-ray system		
Classification panel	Radiology		
Regulation	21 CFR Part 892.1750		
Product Code(s)	JAK		
Legally marketed device(s) to which equivalence is claimed	superDimension i Logic inReach System, K071473/K092365/K102604		
Reason for 510(k) submission	Introduction of accessory Marker Delivery Kit for the placement of markers into soft lung tissue using the superDimension i.Logic inReach System.		
Device description	The superDimension i Logic inReach System is a device that guides a bronchoscope and endoscopic tools to a target in or adjacent to the bronchial tree on a path identified by CT scan. The superDimension i Logic inReach System also allows visualization of the target and the interior of the bronchial tree; placement of catheters in the bronchial tree; and placement of radiosurgical and dye markers into soft lung tissue to guide radiosurgery and thoracic surgery.		
ntended use of the device	Enable the diagnosis and treatment of lung cancer.		

Page 1 of 3

.

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

. •

•

.

dications for use	Indicated for displaying images of the tracheobronchial tree to aid the physician in guiding endoscopic tools or catheters in the pulmonary tract and to enable marker placement within soft lung tissue. It does not make a diagnosis and is not an endoscopic tool. Not for pediatric use.
Summary of the technologic	al characteristics of the device compared to the predicate device
Characteristic	Comparison to New Device
Indications for Use	Same
Anatomical Sites	Same
Basic Principle	Same
Localization Methodology	Same
Visualization Principle and Method of Visualization of Path and Organs	Same
Display Methods	Same
Interventional Instrument	Same. Accessory Marker Delivery Kit is being added.
Method of tracking location	Same
Computer function and type	Same
Energy Released into Body (Negligible)	Same
Compatibility with other devices	Same. Accessory Marker Delivery Kit is being added.
.ocompatibility	Same. Accessory Marker Delivery Kit is being added.
Medical electrical equipment safety compliance	Same
Registration Modes and Navigation Guidance	Same
Software	Same
Sterile Accessories	Same. Accessory Marker Delivery Kit is being added.
Marker Placement	Same. Accessory Marker Delivery Kit is being added.

.

· •

SUMMARY OF NON-CLINICAL TESTS CONDUCTED FOR DETERMINATION OF SUBSTANTIAL EQUIVALENCE

Performance Test Summary-New Device

The superDimension Marker Delivery Kit and Instructions for Use were subjected to the superDimension design control process. Risk Management was performed to analyze the potential hazards associated with the changes. Appropriate design verification and validations were performed to assure the superDimension i Logic inReach System continues to be safe and effective for its intended use.

PERFORMANCE DATA

Comparative Performance Information Summary

Not Required to validate the changes to the superDimension System.

SUMMARY OF CLINICAL TESTS CONDUCTED FOR DETERMINATION OF SUBSTANTIAL EQUIVALENCE AND/OR OF CLINICAL INFORMATION

Clinical tests were not required to validate the changes to the superDimension System.

CONCLUSIONS DRAWN FROM NON-CLINICAL AND CLINICAL DATA

The superDimension i-Logic inReach System with the accessory Marker Delivery Kit is safe and effective for its intended use.

Page 3 of 3



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration 10903 New Hampshire Avenue Document Control Room – WO66-G609 Silver Spring, MD 20993-0002

superDimension, Ltd. % Mr Jonathan Kovach Vice President, Quality and Regulatory Affairs superDimension, Inc. 161 Cheshire Lane, Suite 100 MINNEAPOLIS MN 55441

FEB 1 1 2011

Re: K110093

Trade/Device Name: superDimension® Marker Delivery Kit Regulation Number: 21 CFR 892.1750 Regulation Name: Computed tomography x-ray system Regulatory Class: II Product Code: JAK Dated: January 10, 2011 Received: January 12, 2011

Dear Mr. Kovach:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into class II (Special Controls), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of

Page 2

medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Parts 801 and 809), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm</u> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address http://www.fda.gov/cdrh/industry/support/index.html.

Sincerely Yours,

Mary SVostel

Mary Pastel, ScD. Director Division of Radiological Devices Office of In Vitro Diagnostic Device Evaluation and Safety Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): |< 1/0093

Device Name: superDimension[®] Marker Delivery Kit

Indications for Use:

Indicated for displaying images of the tracheobronchial tree to aid the physician in guiding endoscopic tools or catheters in the pulmonary tract and to enable marker placement within soft lung tissue. It does not make a diagnosis and is not an endoscopic tool. Not for pediatric use.

Prescription Use _____ (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use ______(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

Page 1 of

(Division Sign-Off) Division of Radiological Devices Office of In Vitro Diagnostic Device Evaluation and Safety

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

510(k) Summary superDimension, Ltd. Special 510(k) superDimension® i·Logic[™] inReach® System Addition of Edge[™] Catheter System

OCT - 7 2010 K102604

Date Prepared:

09/07/2010

510(k) Applicant:

superDimension, Ltd. 8 Hamenofim St., P.O. Box 2045 Herzliya 46120 Israel Ph: +972-(0)9-971-3700 Fax: +972-(0)9-971-3701

510(k) Application Correspondent:

Jonathan Kovach Vice President, Quality and Regulatory Affairs superDimension, Inc. 161 Cheshire Lane, Suite 100 Minneapolis, MN 55441 Phone: 763-210-4015 Cell : 763-360-4984 Fax : 763-210-4098 Email : jkovach@superdimension.com

Name of Device :

Trade Name :

superDimension i Logic inReach System superDimension inReach System superDimension/Bronchus

Common Name: Bronchoscope Classification Name: Computed tomography x-ray system 21 CFR Part 892.1750 Product code JAK

Equivalent Legally-Marketed Device:

superDimension i Logic inReach System, K071473/K092365

Description:

The superDimension i Logic inReach System is a device that guides a bronchoscope and endoscopic tools to a target in or adjacent to the bronchial tree on a path identified by CT scan. The superDimension i Logic inReach System also allows visualization of the target and the interior of the bronchial tree; placement of catheters in the bronchial tree; and placement of radiosurgical and dye markers into soft lung tissue to guide radiosurgery and thoracic surgery.

superDimension is introducing the Edge Catheter System for use with the superDimension i-Logic inReach System. The i-Logic inReach System accomodates both the Edge Catheter and the currently available inReach Catheter System.

Intended Use:

Indicated for displaying images of the tracheobronchial tree to aid the physician in guiding endoscopic tools or catheters in the pulmonary tract and to enable marker placement within soft lung tissue. It does not make a diagnosis and is not an endoscopic tool. Not for pediatric use.

Summary of Characteristics Compared to Predicate Device:

The Edge Catheter System is an alternate catheter system for use with the superDimension i Logic inReach System. The Edge Catheter System includes modifications to the existing inReach Catheter System, procedure software, and instructions for use. The Edge Catheter and inReach Catheter systems may both be used with the superDimension i Logic inReach System. No changes are being made to the electromagnetic components or fundamental scientific technology of the i Logic inReach System.

Performance Data:

The Edge Catheter System, i Logic inReach Software, and Instructions for Use were subjected to the superDimension design control process. Risk Management was performed to analyze the potential hazards associated with the changes. Appropriate design verification and validations were performed to assure the superDimension i Logic inReach System continues to be safe and effective for its intended use.

Clinical Data:

Clinical tests were not required to validate the changes to the superDimension inReach System.

Conclusion:

The superDimension i Logic inReach System with the Edge Catheter is safe and effective for its intended use.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration 10903 New Hampshire Avenue Document Control Room – WO66-G609 Silver Spring, MD 20993-0002

SuperDimension, Ltd. % Mr. Jonathan Kovach Vice President, Quality and Regulatory Affairs 161 Cheshire Lane, Suite 100 MINNEAPOLIS MN 55441

OCT - 7 2010

Re: K102604

Trade/Device Name: superDimension[®] i-Logic[™] inReach[®] System with Edge[™] Catheter Regulation Number: 21 CFR 892.1750 Regulation Name: Computed tomography x-ray system Regulatory Class: II Product Code: JAK Dated: September 7, 2010 Received: September 10, 2010

Dear Mr. Kovack:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into class II (Special Controls), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of

medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Parts 801 and 809), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm for the CDRH's Office

of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <u>http://www.fda.gov/cdrh/industry/support/index.html</u>.

Sincerely yours,

David G. Brown, Ph.D. Acting Director Division of Radiological Devices Office of *In Vitro* Diagnostic Device Evaluation and Safety Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known):

K102604

Device Name: superDimension[®] i Logic[™] inReach[®] System with Edge[™] Catheter

Indications for Use:

Indicated for displaying images of the tracheobronchial tree to aid the physician in guiding endoscopic tools or catheters in the pulmonary tract and to enable marker placement within soft lung tissue. It does not make a diagnosis and is not an endoscopic tool. Not for pediatric use.

Prescription Use (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use ______(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)


Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

(Division Sign-Off) Division of Radiological Devices

Office of In Vitro Diagnostic Device Evaluation and Safety

Page 1 of

OCT - 7 2010



Attachment 1: Updated Instructions for Use LungPoint Dilation Balloon

Broncus Medical Inc.

LungPoint[™] Dilation Balloon

REF

Catalog number 10008

Instructions For Use (IFU)

Table of Contents

- 1 Device Description
- 2 Intended Use
- 3 Contraindications
- 4 Warnings / Precautions
- 5 Possible Complications
- 6 Preparation for Use
- 7 Operating Instructions
- 8 Storage

Broncus Medical Inc.

RonlyCaution: Federal law (USA) restricts this device to by or on the order of a physician.	
STERILE R	Contents are supplied irradiation sterilized.
2	For single-patient-use only. Do not re-use or re-sterilize.

1 Device Description

The LungPoint Dilation Balloon is an endoscopic tool that is used to dilate the target lung tissue of the bronchial tree. The dilation balloon is inserted through the bronchoscope or other endoscopic tool, such as a sheath, and is used to dilate lung tissue. The balloon may be used during interventional bronchoscopy procedures as determined by the physician performing the bronchoscopy.

The LungPoint Dilation Balloon is comprised of a catheter with a non-compliant 4mm OD x 6mm long inflatable balloon at the distal end and a y-connector with female luer fittings at the proximal end. The side port of the Y connector is used for inflation and deflation of the balloon with a standard indeflator capable of 12atm pressure (typically with a gage) and a male luer lock connector. The center port contains a nitinol wire which is fixed in place for pushability and kink resistance during tortuous bronchoscopic navigation. Included are two radiopaque marker bands and the distal and proximal ends of the balloon inflation length for visualization of the device position under fluoroscopy

Product Configuration and Dimensions

The LungPoint Dilation Balloon (see illustration in Figure 1) has the following specifications:

Catalog	Balloon Size	Rated Balloon	Maximum Catheter	Catheter
Number	(OD x length)	Pressure	Outer Diameter (OD)	Length/working length
10008-1	4mm x 6mm	<u>20atm</u>	1mm	1430mm/975mm

Broncus Medical Inc.

2 Intended Use

The LungPoint Dilation Balloon is part of the LungPoint Tools, which are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs in a minimally invasive manner.

3 Contraindications

There are no known contraindications.




Figure 1: Dilation Balloon (shown inflated)

(I) WARNINGS AND PRECAUTIONS

4 Warnings / Precautions

- The LungPoint Dilation Balloon should only be used by physicians thoroughly trained in bronchoscopy and in the particular technique and procedure to be performed, and familiar with the associated risks. These operating instructions must be read, understood, and followed.
- The LungPoint Dilation Balloon should only be used in a manner consistent with the Bronchoscope manufacturer's Instructions for Use. Use of the LungPoint Dilation Balloon in contradiction to the instructions may harm the device, bronchoscope or patient.
- The rated burst pressure is 20 atmosphere, exceeding the rated burst pressure may result in harm to the patient.

B

Broncus Medical Inc.

5 Possible Complications

- Bleeding
- Infection

6 Preparation for Use

- 6.1 Materials Required
 - LungPoint Dilation Balloon
 - Indeflator capable of 12atm
 - Flexible bronchoscope with a working channel of 2.0 mm or greater
- 6.2 Device Inspection and Preparation
 - 6.2.1 Read all IFUs, package inserts, labels, and warnings for the LungPoint Dilation Balloon before beginning the procedure.

CAUTION: Inspect the sealed sterile device package before opening. If the seal is broken, contents may not be sterile and could pose a risk of patient infection.

6.2.2 Aseptically remove the LungPoint Dilation Balloon from the package and inspect for any damage, such as broken or crushed areas of the catheter shaft, sharp or protruding edges at the distal tip, or any kinks in the catheter shaft.

CAUTION: If any damage or irregularity is found, **DO NOT USE** the device as this could result in harm to the patient.

CAUTION: DO NOT USE the device past its "Use Before" (expiration) date.

CAUTION: DO NOT place or manipulate the catheter in the airways unless under direct visualization with a bronchoscope or using fluoroscopy as this may result in harm to the patient.

Broncus Medical Inc.

7 Operating Instructions

7.1 Carefully advance the balloon through the working channel of the bronchoscope using short (approximately 20 mm) increments, until the distal end of the catheter is within bronchoscopic view.

CAUTION: DO NOT force the catheter if resistance to insertion is encountered. Reduce the angulation of the bronchoscope until the device passes smoothly. This can result in (a) kinks to the device, (b) damage to the bronchoscope, or (c) harm to the patient such as punctures, hemorrhage or mucous membrane damage.

CAUTION: If any damage occurs, including a kink, DO NOT USE and discard.

7.2 Advance the catheter to the desired protrusion from the bronchoscope as required for navigation to the region of interest.

CAUTION: Do not angulate the bronchoscope abruptly while the device is extended from the distal end of the bronchoscope. This can result in kinks to the device or harm to the patient such as punctures, hemorrhage or mucous membrane damage.

- 7.3 Connect the indeflator via the luer lock connector and inflate to 10atm. Do not exceed the rated burst pressure of 20 atm.
- 7.4 If using fluoroscopy visualize balloon inflation location via the radiopaque marker bands and the distal and proximal ends of the balloon inflation length.
- 7.5 Deflate the balloon, and remove the device

CAUTION: If excessive resistance makes withdrawal difficult, reduce the angulation of the bronchoscope until the accessory can be withdrawn smoothly. Forcible withdrawal could damage the device or bronchoscope.

7.6 At the completion of the patient procedure, dispose of the device in accordance with applicable hospital, local, state and federal laws and regulations.

NOTE: The LungPoint Dilation Balloon has been tested to three insertion cycles, meaning that the dilation balloon can be inserted/inflated/deflated/removed a total of three times and must be discarded after the third complete cycle.

8 Storage

Store at controlled room temperature.

😽 Broncus Medical Inc.

LEGAL NOTICE

© 2013 Broncus Medical, Inc. All rights reserved.

Broncus, the Broncus logo, and LungPoint are trademarks or registered trademarks of Broncus Medical, Inc.

Patents Pending.

REF	Catalog number	
Ronly	Federal law (USA) restricts this device to sale by or on the order of a physician.	
<u></u>	Caution: consult instructions for use	
STERILE R	Irradiation sterilized. Sterility guaranteed if package unopened and undamaged.	
LOT Serial number or batch code		
Use before		
2 For single use		
	Manufacturer	

Graphic Symbol Legend for Medical Device Labeling



Broncus Medical Inc.

Manufacturer: Broncus Medical, Inc. 1400 N. Shoreline Blvd., Bldg. A, Suite 8 Mountain View, CA, 94043 USA (877) 428-1600

Attachment 2: Updated 510(k) Summary

510(k) Summary

1. Date of Summary

April 22, 2013

2. 510(k) Applicant

Broncus Medical, Inc. 1400 N. Shoreline Blvd, Suite A8 Mountain View, California 94043 Phone: (650) 428-1600 FAX: (650) 428-1542

Contact Person:	Gary Kaplan
Phone:	(650) 428-1600
Fax:	(650) 428-1542
e-mail:	gkaplan@broncus.com

3. Device Overview

Trade Name:	LungPoint [™] Tools (LungPoint Sheath and LungPoint Dilation Balloon)
Common Name:	Sheath and Dilation Balloon
Classification Name:	Bronchoscope and Accessories 21 CFR 874.4680 Product Code EOQ

4. Predicate Device

The predicate devices identified are as follows:

Trade Name	510(k) Submitter	510(k) Number		
LungPoint Planning and Virtual Bronchoscopic Navigation (VBN) System (with FlexNeedle)	Broncus Technologies, Inc*	K112051, cleared to market on October 12, 2011 (and K093423, K091160 and K090095, by reference)		
inReach System	SuperDimension	K110093, cleared to market on February 11, 2011 (and K071473, K092365 and K102604, by reference)		
CRE Pulmonary Balloon Dilatation Catheter	Boston Scientific CorporationK023337, cleared to marke November 18, 2002			
*NOTE: Broncus Technologies, Inc changed their name to Broncus Medical, Inc in June 2012.				

Broncus Medical, Inc.

Page 1 of 3

5. Device Description

The LungPoint Tools are endoscopic tools used during bronchoscopy procedures. The LungPoint Sheath is designed to be used with a bronchoscope to provide a working channel through which endoscopic tools, such as needles, dilation balloons, or other endoscopic devices may be introduced to the targeted lung tissue within the respiratory organs. It is advanced to a predefined target following guidance of the LungPoint Software. The LungPoint Dilation Balloon is used to dilate tissue of the bronchial tree and may be inserted through the sheath or directly through the working channel of the bronchoscope.

The materials used in the LungPoint Tools are commonly used medical grade materials and include platinum iridium markers.

Catalog Number	Working Length	Maximum Catheter Outer Diameter (OD)	Catheter Internal Diameter (ID)	Minimum accessory length for use through sheath
10007-1	900mm	2.65mm	2.0mm	965mm

The LungPoint Sheath has the following specifications:

The LungPoint Dilation Balloon has the following specifications:

Catalog Number	Balloon Size (OD x length)	Rated Balloon Pressure	Maximum Catheter Outer Diameter (OD)	Catheter Length/working length
10008-1	4mm x 6mm	20atm	1mm	1430mm/975mm

6. Intended Use

The LungPoint Tools are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs in a minimally invasive manner.

7. Comparison to Predicate Device

The LungPoint Tools are commonly used endoscopic tools with the same technological characteristics as the predicate devices. The indications for use of the LungPoint Tools for use with the LungPoint Software are all within the intended use of the predicate device, which is to enable treatment and diagnosis of lung cancer. The technological

characteristics of the subject devices are the same as those of the predicate devices with the exception of the balloon size and length, balloon burst pressure and balloon catheter length as outlined.

	CRE Pulmonary Balloon Dilatation Catheter (K023337)	LungPoint Dilation Balloon
Balloon Size and Inflation	8mm @ 3 ATM	4mm @ 10 ATM
(balloon pressure)	9mm @ 5.5 ATM	
	10mm @ 9 ATM	
Catheter Length (cm)	155	143
Balloon Length (cm)	3.0	0.6
Rated Burst Pressure (atmospheres)	9	20

None of these differences raise new questions of safety and effectiveness. Performance of the subject devices has been verified by use of accepted methods.

8. Performance Data

The design and safety of the LungPoint Tools were verified by performing functional and performance testing. All tests were designed to subject the sheath and dilation balloon to stresses that exceed those which would be encountered during clinical use. Testing included the following:

- Dimensional testing
- Joint/tensile test
- Simulated use
- Balloon fatigue/burst pressure
- Radiopacity.

All testing results met the pre-determined acceptance criteria that were established in the test protocols. Based on the testing the LungPoint Tools are as safe, as effective, and perform at least as safely and effectively as the predicate devices.

9. Safety and Effectiveness

The LungPoint Tools labeling contains instructions for use and any necessary cautions and warnings, to provide for safe and effective use of the devices. The biocompatibility assessment of all patient contacting materials was performed in accordance with ISO 10993, *Biological Evaluation of Medical Devices*. Specifically, cytotoxicity, sensitization, intracutaneous reactivity and systemic toxicity (acute) were tested. In addition, the devices are sterilized using e-beam sterilization.

10. Conclusion

Based on the testing the LungPoint Tools are as safe, as effective, and perform at least as safely and effectively as the predicate devices.

Attachment 3: FDA Form 3654 for ISO 15223

Records processed under FOIA Request # 2015-7253	<u>}.FREI&&BEUBY@DRH9 88102169-2016</u>	ation Dat	te: 12/31/13			
Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)						
This report and the Summary Report Table are to be completed by ences a national or international standard. A separate report is requi						
TYPE OF 510(K) SUBMISSION	bbreviated					
STANDARD TITLE ¹ BS EN ISO 15223-1:2012 Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied						
Please answer the following questions		Yes	No			
Is this standard recognized by FDA ² ?		\times				
FDA Recognition number ³	# <u>5</u>	5-73				
Was a third party laboratory responsible for testing conformity of the in the 510(k)?			\times			
Is a summary report ⁴ describing the extent of conformance of the standard used included in the 510(k)? If no, complete a summary report table.			X			
Does the test data for this device demonstrate conformity to the requirements to this device?	X					
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).			X			
Does this standard include more than one option or selection of test If yes, report options selected in the summary report table.	s?	\boxtimes				
Were there any deviations or adaptations made in the use of the sta If yes, were deviations in accordance with the FDA supplemental inf						
Were deviations or adaptations made beyond what is specified in the lf yes, report these deviations or adaptations in the summary report			\boxtimes			
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			\times			
Is there an FDA guidance ⁶ that is associated with this standard? If yes, was the guidance document followed in preparation of this 51 Title of guidance:						
[title of standard] [date of publication] standa 2 Authority [21 U.S.C. 360d], www.fda.gov/cdrh/stdsprog.html utilized 3 http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/ 5 The su which i search.cfm 4 The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); 6 The on choices made when options or a selection of methods are descr bed; 6 The on	ation body involved in conformance assessment and. The summary report includes information on a d during the development of the device. upplemental information sheet (SIS) is additional i is necessary before FDA recognizes the standard www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfSta n.cfm nline search for CDRH Guidance Documents can da.gov/cdrh/guidance.html	all standa informati d. Founc andards/	on I at			

FORM FDA 3654e(12/10)? Contact FDA/CDRH/OCE/DID at PORH-FOISTATUS@fda.hhs.gov or 301-796-8198 Phics (301) 443-6740 EF

EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE

SUMMART REPORT TABLE					
STANDARD TITLE BS EN ISO 15223-1:2012 Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied					
CONFORMANCE WITH STANDARD SECTIONS*					
SECTION NUMBER	SECTION TITLE		CONFORMANCE?		
1, 2, 3	Scope, Normative references, Terms and defin	itions	🗙 Yes 🗌 No 🗌 N/A		
TYPE OF DEVIATION OF NA	ROPTION SELECTED *				
DESCRIPTION NA					
JUSTIFICATION NA					
SECTION NUMBER	SECTION TITLE		CONFORMANCE?		
4	General Requirements		🗙 Yes 🗌 No 🗌 N/A		
TYPE OF DEVIATION OR OPTION SELECTED * NA					
DESCRIPTION NA					
JUSTIFICATION NA					
SECTION NUMBER	SECTION TITLE		CONFORMANCE?		
5	Symbols		🗙 Yes 🗌 No 🗌 N/A		
TYPE OF DEVIATION OR OPTION SELECTED * Symbols in the following Sections were applied: 5.1.1, 5.1.4, 5.1.5, 5.1.6, 5.1.7, 5.2.4, 5.4.2, 5.4.3					
	rom the sections above.				
JUSTIFICATION Not all symbols in the standard were applicable.					
* For completeness list all sections of the standard and indicate whether conformance is met. If a section is not applicable (N/A) an explanation is needed under "justification." Some standards include options, so similar to deviations, the option chosen needs to be described and adequately justified as appropriate for the subject device. Explanation of all deviations or description of options selected when following a standard is required under "type of deviation or option selected," "description" and "justification" on the report. More than one page may be necessary.					
* Types of deviations can include an exclusion of a section in the standard, a deviation brought out by the FDA supplemental information sheet (SIS), a deviation to adapt the standard to the device, or any adaptation of a section.					
Paperwork Reduction Act Statement					
Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to:					
Food a Office 1350 P	ment of Health and Human Services nd Drug Administration of Chief Information Officer ficcard Drive, Room 400 ille, MD 20850	An agency may not conduct or spon required to respond to, a collection displays a currently valid OMB con	of information unless it		



Broncus Medical, Inc. (650) 428-1600 1400 N. Shoreline Blvd. (650) 428-1542 fax Suite A8 Mountain View, CA 94043 www.broncus.com

April 22, 2013

Dear Ms. (b)(6)

In response to your email dated April 19, 2013, we are submitting the attached response to your questions. Your questions are identified in **bold** text. The corresponding Broncus responses are identified in normal text immediately following each question or part of a question.

The following attachments are included with the response,

Attachment 1: Updated Instructions for Use LungPoint Dilation Balloon Attachment 2: Updated 510(k) Summary Attachment 3: FDA Form 3654 for ISO 15223.

With the exception of the 510(k) summary, changes to the attached files, where appropriate, are underlined.

Please feel free to contact me should you have any additional questions.

Regards,



Attachment 1: Updated Instructions for Use LungPoint Dilation Balloon

Broncus Medical Inc.

LungPoint[™] Dilation Balloon

REF

Catalog number 10008

Instructions For Use (IFU)

Table of Contents

- 1 Device Description
- 2 Intended Use
- 3 Contraindications
- 4 Warnings / Precautions
- 5 Possible Complications
- 6 Preparation for Use
- 7 Operating Instructions
- 8 Storage

Broncus Medical Inc.

Ronly	Caution : Federal law (USA) restricts this device to sale by or on the order of a physician.
STERILE R	Contents are supplied irradiation sterilized.
2	For single-patient-use only. Do not re-use or re-sterilize.

1 Device Description

The LungPoint Dilation Balloon is an endoscopic tool that is used to dilate the target lung tissue of the bronchial tree. The dilation balloon is inserted through the bronchoscope or other endoscopic tool, such as a sheath, and is used to dilate lung tissue. The balloon may be used during interventional bronchoscopy procedures as determined by the physician performing the bronchoscopy.

The LungPoint Dilation Balloon is comprised of a catheter with a non-compliant 4mm OD x 6mm long inflatable balloon at the distal end and a y-connector with female luer fittings at the proximal end. The side port of the Y connector is used for inflation and deflation of the balloon with a standard indeflator capable of 12atm pressure (typically with a gage) and a male luer lock connector. The center port contains a nitinol wire which is fixed in place for pushability and kink resistance during tortuous bronchoscopic navigation. Included are two radiopaque marker bands and the distal and proximal ends of the balloon inflation length for visualization of the device position under fluoroscopy

Product Configuration and Dimensions

The LungPoint Dilation Balloon (see illustration in Figure 1) has the following specifications:

Catalog	Balloon Size	Rated Balloon	Maximum Catheter	Catheter
Number	(OD x length)	Pressure	Outer Diameter (OD)	Length/working length
10008-1	4mm x 6mm	<u>20atm</u>	1mm	1430mm/975mm

Broncus Medical Inc.

2 Intended Use

The LungPoint Dilation Balloon is part of the LungPoint Tools, which are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs in a minimally invasive manner.

3 Contraindications

There are no known contraindications.




Figure 1: Dilation Balloon (shown inflated)

(I) WARNINGS AND PRECAUTIONS

4 Warnings / Precautions

- The LungPoint Dilation Balloon should only be used by physicians thoroughly trained in bronchoscopy and in the particular technique and procedure to be performed, and familiar with the associated risks. These operating instructions must be read, understood, and followed.
- The LungPoint Dilation Balloon should only be used in a manner consistent with the Bronchoscope manufacturer's Instructions for Use. Use of the LungPoint Dilation Balloon in contradiction to the instructions may harm the device, bronchoscope or patient.
- The rated burst pressure is 20 atmosphere, exceeding the rated burst pressure may result in harm to the patient.

B

Broncus Medical Inc.

5 Possible Complications

- Bleeding
- Infection

6 Preparation for Use

- 6.1 Materials Required
 - LungPoint Dilation Balloon
 - Indeflator capable of 12atm
 - Flexible bronchoscope with a working channel of 2.0 mm or greater
- 6.2 Device Inspection and Preparation
 - 6.2.1 Read all IFUs, package inserts, labels, and warnings for the LungPoint Dilation Balloon before beginning the procedure.

CAUTION: Inspect the sealed sterile device package before opening. If the seal is broken, contents may not be sterile and could pose a risk of patient infection.

6.2.2 Aseptically remove the LungPoint Dilation Balloon from the package and inspect for any damage, such as broken or crushed areas of the catheter shaft, sharp or protruding edges at the distal tip, or any kinks in the catheter shaft.

CAUTION: If any damage or irregularity is found, **DO NOT USE** the device as this could result in harm to the patient.

CAUTION: DO NOT USE the device past its "Use Before" (expiration) date.

CAUTION: DO NOT place or manipulate the catheter in the airways unless under direct visualization with a bronchoscope or using fluoroscopy as this may result in harm to the patient.

Broncus Medical Inc.

7 Operating Instructions

7.1 Carefully advance the balloon through the working channel of the bronchoscope using short (approximately 20 mm) increments, until the distal end of the catheter is within bronchoscopic view.

CAUTION: DO NOT force the catheter if resistance to insertion is encountered. Reduce the angulation of the bronchoscope until the device passes smoothly. This can result in (a) kinks to the device, (b) damage to the bronchoscope, or (c) harm to the patient such as punctures, hemorrhage or mucous membrane damage.

CAUTION: If any damage occurs, including a kink, DO NOT USE and discard.

7.2 Advance the catheter to the desired protrusion from the bronchoscope as required for navigation to the region of interest.

CAUTION: Do not angulate the bronchoscope abruptly while the device is extended from the distal end of the bronchoscope. This can result in kinks to the device or harm to the patient such as punctures, hemorrhage or mucous membrane damage.

- 7.3 Connect the indeflator via the luer lock connector and inflate to 10atm. Do not exceed the rated burst pressure of 20 atm.
- 7.4 If using fluoroscopy visualize balloon inflation location via the radiopaque marker bands and the distal and proximal ends of the balloon inflation length.
- 7.5 Deflate the balloon, and remove the device

CAUTION: If excessive resistance makes withdrawal difficult, reduce the angulation of the bronchoscope until the accessory can be withdrawn smoothly. Forcible withdrawal could damage the device or bronchoscope.

7.6 At the completion of the patient procedure, dispose of the device in accordance with applicable hospital, local, state and federal laws and regulations.

NOTE: The LungPoint Dilation Balloon has been tested to three insertion cycles, meaning that the dilation balloon can be inserted/inflated/deflated/removed a total of three times and must be discarded after the third complete cycle.

8 Storage

Store at controlled room temperature.

😽 Broncus Medical Inc.

LEGAL NOTICE

© 2013 Broncus Medical, Inc. All rights reserved.

Broncus, the Broncus logo, and LungPoint are trademarks or registered trademarks of Broncus Medical, Inc.

Patents Pending.

REF	Catalog number
Ronly	Federal law (USA) restricts this device to sale by or on the order of a physician.
<u></u>	Caution: consult instructions for use
STERILE R	Irradiation sterilized. Sterility guaranteed if package unopened and undamaged.
LOT	Serial number or batch code
	Use before
2	For single use
	Manufacturer

Graphic Symbol Legend for Medical Device Labeling



Broncus Medical Inc.

Manufacturer: Broncus Medical, Inc. 1400 N. Shoreline Blvd., Bldg. A, Suite 8 Mountain View, CA, 94043 USA (877) 428-1600

Attachment 2: Updated 510(k) Summary

510(k) Summary

1. Date of Summary

April 22, 2013

2. 510(k) Applicant

Broncus Medical, Inc. 1400 N. Shoreline Blvd, Suite A8 Mountain View, California 94043 Phone: (650) 428-1600 FAX: (650) 428-1542

Contact Person:	Gary Kaplan
Phone:	(650) 428-1600
Fax:	(650) 428-1542
e-mail:	gkaplan@broncus.com

3. Device Overview

Trade Name:	LungPoint [™] Tools (LungPoint Sheath and LungPoint Dilation Balloon)
Common Name:	Sheath and Dilation Balloon
Classification Name:	Bronchoscope and Accessories 21 CFR 874.4680 Product Code EOQ

4. Predicate Device

The predicate devices identified are as follows:

Trade Name	510(k) Submitter	510(k) Number		
LungPoint Planning and Virtual Bronchoscopic Navigation (VBN) System (with FlexNeedle)	Broncus Technologies, Inc*	K112051, cleared to market on October 12, 2011 (and K093423, K091160 and K090095, by reference)		
inReach System	SuperDimension	K110093, cleared to market on February 11, 2011 (and K071473, K092365 and K102604, by reference)		
CRE Pulmonary Balloon Dilatation CatheterBoston Scientific CorporationK023337, cleared to market on November 18, 2002				
*NOTE: Broncus Technologies, Inc ch	anged their name to Bron	cus Medical, Inc in June 2012.		

5. Device Description

The LungPoint Tools are endoscopic tools used during bronchoscopy procedures. The LungPoint Sheath is designed to be used with a bronchoscope to provide a working channel through which endoscopic tools, such as needles, dilation balloons, or other endoscopic devices may be introduced to the targeted lung tissue within the respiratory organs. It is advanced to a predefined target following guidance of the LungPoint Software. The LungPoint Dilation Balloon is used to dilate tissue of the bronchial tree and may be inserted through the sheath or directly through the working channel of the bronchoscope.

The materials used in the LungPoint Tools are commonly used medical grade materials and include platinum iridium markers.

Catalog Number	Working Length	Maximum Catheter Outer Diameter (OD)	Catheter Internal Diameter (ID)	Minimum accessory length for use through sheath
10007-1	900mm	2.65mm	2.0mm	965mm

The LungPoint Sheath has the following specifications:

The LungPoint Dilation Balloon has the following specifications:

Catalog Number	Balloon Size (OD x length)	Rated Balloon Pressure	Maximum Catheter Outer Diameter (OD)	Catheter Length/working length
10008-1	4mm x 6mm	20atm	1mm	1430mm/975mm

6. Intended Use

The LungPoint Tools are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs in a minimally invasive manner.

7. Comparison to Predicate Device

The LungPoint Tools are commonly used endoscopic tools with the same technological characteristics as the predicate devices. The indications for use of the LungPoint Tools for use with the LungPoint Software are all within the intended use of the predicate device, which is to enable treatment and diagnosis of lung cancer. The technological

Appendix A

characteristics of the subject devices are the same as those of the predicate devices with the exception of the balloon size and length, balloon burst pressure and balloon catheter length as outlined.

	CRE Pulmonary Balloon Dilatation Catheter (K023337)	LungPoint Dilation Balloon
Balloon Size and Inflation	8mm @ 3 ATM	4mm @ 10 ATM
(balloon pressure)	9mm @ 5.5 ATM	
	10mm @ 9 ATM	
Catheter Length (cm)	155	143
Balloon Length (cm)	3.0	0.6
Rated Burst Pressure (atmospheres)	9	20

None of these differences raise new questions of safety and effectiveness. Performance of the subject devices has been verified by use of accepted methods.

8. Performance Data

The design and safety of the LungPoint Tools were verified by performing functional and performance testing. All tests were designed to subject the sheath and dilation balloon to stresses that exceed those which would be encountered during clinical use. Testing included the following:

- Dimensional testing
- Joint/tensile test
- Simulated use
- Balloon fatigue/burst pressure
- Radiopacity.

All testing results met the pre-determined acceptance criteria that were established in the test protocols. Based on the testing the LungPoint Tools are as safe, as effective, and perform at least as safely and effectively as the predicate devices.

9. Safety and Effectiveness

The LungPoint Tools labeling contains instructions for use and any necessary cautions and warnings, to provide for safe and effective use of the devices. The biocompatibility assessment of all patient contacting materials was performed in accordance with ISO 10993, *Biological Evaluation of Medical Devices*. Specifically, cytotoxicity, sensitization, intracutaneous reactivity and systemic toxicity (acute) were tested. In addition, the devices are sterilized using e-beam sterilization.

10. Conclusion

Based on the testing the LungPoint Tools are as safe, as effective, and perform at least as safely and effectively as the predicate devices.

Appendix A Page 3 of 3

Attachment 3: FDA Form 3654 for ISO 15223

Records processed under FOIA Request # 2015-725	3.FREI&8880195/2015RH9 88102169-2091	gtion Da	te: 12/31/13
Department of Health and Hu Food and Drug Admini STANDARDS DATA REPO (To be filled in by ap	istration IRT FOR 510(k)s		
This report and the Summary Report Table are to be completed be ences a national or international standard. A separate report is requ			
TYPE OF 510(K) SUBMISSION	Abbreviated		
STANDARD TITLE ¹ BS EN ISO 15223-1:2012 Medical devices - Symbols to be used with med	ical device labels, labeling and informati	ion to be	e supplied
Please answer the following questions		Yes	No
Is this standard recognized by FDA ² ?		\times	
FDA Recognition number ³		5-73	
Was a third party laboratory responsible for testing conformity of th in the 510(k)?			\times
Is a summary report ⁴ describing the extent of conformance of the s 510(k)? If no, complete a summary report table.			\boxtimes
Does the test data for this device demonstrate conformity to the rec pertains to this device?		\boxtimes	
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).			\times
Does this standard include more than one option or selection of tes If yes, report options selected in the summary report table.	sts?	\times	
Were there any deviations or adaptations made in the use of the st If yes, were deviations in accordance with the FDA supplemental ir			\boxtimes
Were deviations or adaptations made beyond what is specified in the lf yes, report these deviations or adaptations in the summary report			\boxtimes
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			\boxtimes
Is there an FDA guidance ⁶ that is associated with this standard? If yes, was the guidance document followed in preparation of this 5 Title of guidance:			
[title of standard] [date of publication] stand 2 Authority [21 U.S.C. 360d], www.fda.gov/cdrh/stdsprog.html utilize 3 http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/ 5 The s search.cfm which 4 The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); search secret choices made when options or a selection of methods are descr bed; 6 The comparison	ication body involved in conformance assessment lard. The summary report includes information on ed during the development of the device. supplemental information sheet (SIS) is additional in is necessary before FDA recognizes the standar /www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfSi sh.cfm online search for CDRH Guidance Documents car fda.gov/cdrh/guidance.html	all stand informati rd. Found tandards/	on d at

FORM FDA 3654e(12/10)? Contact FDA/CDRH/OCE/DID at PORH-FOISTATUS@fda.hhs.gov or 301-796-8198 Phics (301) 443-6740 EF

EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE

	SOWIWART REF		
STANDARD TITLE BS EN ISO 15223-1:20	012 Medical devices - Symbols to be used with	medical device labels, labeling a	nd information to be supplied
	CONFORMANCE WITH STA	NDARD SECTIONS*	
SECTION NUMBER	SECTION TITLE		CONFORMANCE?
1, 2, 3	Scope, Normative references, Terms and defin	itions	🗙 Yes 🗌 No 🗌 N/A
TYPE OF DEVIATION OF NA	ROPTION SELECTED *		
DESCRIPTION NA			
JUSTIFICATION NA			
SECTION NUMBER	SECTION TITLE		CONFORMANCE?
4	General Requirements		🗙 Yes 🗌 No 🗌 N/A
TYPE OF DEVIATION OF NA	R OPTION SELECTED *		
DESCRIPTION NA			
JUSTIFICATION NA			
SECTION NUMBER	SECTION TITLE		CONFORMANCE?
5	Symbols		🗙 Yes 🗌 No 🗌 N/A
TYPE OF DEVIATION OF Symbols in the followin	R OPTION SELECTED * ng Sections were applied: 5.1.1, 5.1.4, 5.1.5, 5.1	1.6, 5.1.7, 5.2.4, 5.4.2, 5.4.3	
	rom the sections above.		
JUSTIFICATION Not all symbols in the s	standard were applicable.		
explanation is needed described and adequa selected when followi report. More than on	all sections of the standard and indicate whether d under "justification." Some standards include of ately justified as appropriate for the subject devi ng a standard is required under "type of deviation e page may be necessary.	options, so similar to deviations, th ice. Explanation of all deviations o on or option selected," "description	e option chosen needs to be r description of options " and "justification" on the
	an include an exclusion of a section in the stand S), a deviation to adapt the standard to the devi		
	Paperwork Reduction	Act Statement	
time for reviewi completing and	burden for this collection of information is est ng instructions, searching existing data sources reviewing the collection of information. Send collection of information, including suggestions f	s, gathering and maintaining the decomments regarding this burden estimates and the state of the	ata needed, and
Food a Office 1350 P	ment of Health and Human Services nd Drug Administration of Chief Information Officer ficcard Drive, Room 400 ille, MD 20850	An agency may not conduct or spon required to respond to, a collection displays a currently valid OMB con	of information unless it



lana an bhannai, but 1,00 M Storiana Bas Easter Aub Moonta 5 Rev. CA 24123 a de proto a clam

3656 A18 1600 logeration tax K131234/51

August 20, 2013



In response to your email dated June 12, 2013, we are submitting the attached response to the questions posed by FDA. FDA's questions are identified in bold text. The corresponding Broncus responses are identified in normal text immediately following each question.

The following attachments are included with the response,

Attachment 1: Updated Section 12 - Substantial Equivalence Discussion Attachment 2: Updated 510(k) Summary Attachment 3: Predicate Device Labeling - Olympus Guide Sheath Attachment 4: Instructions for Use LungPoint Dilation Balloon Attachment 5: Instructions for Use LungPoint Sheath Attachment 6:(b)(4) Attachment 7: Data Attachment 8: Attachment 9:

This submission includes one original paper copy and one electronic copy. The electronic copy included here is an exact duplicate of the original paper submission.

Please feel	free to	contact n	ne should	you h	ave any	additional	questions at
(b)(4)							

. . .

Regards,



KA CONSULTAIL

Broncus Medical, Inc. CONFIDENTIAL

Page 1 of 12



Broncus Medical, Inc. 1400 N. Shoreline Blvd. Suite A8 Mountain View, CA 94043 www.broncus.com (650) 428-1600 (650) 428-1542 fax

August 20, 2013

Dear Ms (b)(6)

In response to your email dated June 12, 2013, we are submitting the attached response to the questions posed by FDA. FDA's questions are identified in bold text. The corresponding Broncus responses are identified in normal text immediately following each question.

The following attachments are included with the response,

Attachment 1: Updated Section 12 - Substantial Equivalence Discussion Attachment 2: Updated 510(k) Summary Attachment 3: Predicate Device Labeling - Olympus Guide Sheath Attachment 4: Instructions for Use LungPoint Dilation Balloon Attachment 5: Instructions for Use LungPoint Sheath Attachment 6 Attachment 7 Data Attachment 8 Attachment 8 Attachment 9

This submission includes one original paper copy and one electronic copy. The electronic copy included here is an exact duplicate of the original paper submission.

Please feel free to	o contact m	ne should	you have	any additional	questions at
(b)(4)				2	

	Regards,	
(b)(6)		

Attachment 1: Updated Section 12 - Substantial Equivalence Discussion

SECTION 12: SUBSTANTIAL EQUIVALENCE DISCUSSION

12.1 Predicate Devices

The LungPoint Tools (LungPoint Dilation Balloon and LungPoint Sheath), which are designed to be used as accessories to the LungPoint Software System (cleared to market with a transbronchial aspiration needle on October 12, 2011 (K112051)), are commonly used endoscopic tools.

The predicate devices identified for the LungPoint Dilation Balloon and LungPoint Sheath are as follows:

Trade Name	510(k) Submitter	510(k) Number
Olympus Guide Sheath	Olympus Medical Systems Corporation	K060243, cleared to market on June 23, 2006
CRE Pulmonary Balloon Dilatation Catheter	Boston Scientific Corporation	K023337, cleared to market on November 18, 2002

12.2 Comparison to the Predicate Devices

The following tables outline the similarities and differences between the LungPoint Tools and the predicate devices.

14010 12.2. 1. 0	Table 12.21: Comparison Table - LungPoint Sneath		
	Olympus Guide Sheath	LungPoint Sheath	
	(K060243)		
Indications for use	This instrument has been designed to be used with Olympus bronchoscopes, endo-therapy accessories and ultrasound probe to guide the endo-therapy accessories or ultrasound probe to the target area within the respiratory organs.	The LungPoint Tools are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs in a minimally invasive manner.	
(4)			
Fluoro Tip/Radiopaque	Yes	Yes	
Markers			
Delivery Approach	Visual (via bronchoscope) and/or Fluoroscopy	Visual (via bronchoscope) and/or Fluoroscopy	

Yes

Table 12.21: Compariso	n Table - LungPoint Sheath
------------------------	----------------------------

Single Use

Yes

		CRE Pulmonary Balloon Dilatation Catheter	LungPoint Dilation Balloon	
		(K023337)		
	Indications for use	The CRE Pulmonary Balloon Dilatation Catheter is intended to be used endoscopically to dilate strictures of the airway tree.	The LungPoint Tools are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs in a minimally invasive manner.	
	Tensile Strength	Complies with ISO 10555 as it is the industry standard	Complies with ISO 10555	
	Delivery Approach	Visual via bronchoscope	Visual via bronchoscope	
	Delivered Through	Flexible bronchoscope with minimum working channel of mm	Flexible bronchoscope with minimum working channel of b.mm	
	Single Use	Yes	Yes	
(b)(4)				

Table 12.2-2: Comparison Table - LungPoint Dilation Balloon

12.3 Substantial Equivalence Discussion


The LungPoint Tools (LungPoint Sheath and Dilation Balloon), the subject devices of this premarket notification, are substantially equivalent to the predicate devices. As outlined above, the LungPoint Tools have the same intended use and technological characteristics as their respective predicate devices.

In making this determination, Broncus followed the FDA's Draft Guidance entitled The 510(k) Program: Evaluating Substantial Equivalence in Premarket Notifications [510(k)] (December 27, 2011), the detailed flowchart provided in the guidance (included at the end of this section) and 21 CFR 807.100.

Specifically the following questions were considered:


b)(4)

(b)(4) Deficiencies	
(b)(4) Deficiencies	



CONFIDENTIAL Broncus Medical, Inc.

Page 7 of 8 Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 30ATtecRintent 1 Page 8 of 9



This Flowchart is not intended to be used as a 'stand-alone' document and should only be considered in conjunction with the accompanying text in this guidance.

Figure 1: 510(k) "Substantial Equivalence" Decision-Making Flowchart

Attachment 2: Updated 510(k) Summary

510(k) Summary

1. Date of Summary

August 20, 2013

2. 510(k) Applicant

Broncus Medical, Inc. 1400 N. Shoreline Blvd, Suite A8 Mountain View, California 94043 Phone: (650) 428-1600 FAX: (650) 428-1542

Contact Person:	Gary Kaplan
Phone:	(650) 428-1600
Fax:	(650) 428-1542
e-mail:	<u>gkaplan@broncus.com</u>

3. Device Overview

Trade Name:	LungPoint [™] Tools (LungPoint Sheath and LungPoint Dilation Balloon)
Common Name:	Sheath and Dilation Balloon
Classification Name:	Bronchoscope and Accessories 21 CFR 874.4680 Product Code EOQ

4. Predicate Device

The predicate devices identified are as follows:

Trade Name	510(k) Submitter	510(k) Number
CRE Pulmonary Balloon Dilatation Catheter	Boston Scientific Corporation	K023337, cleared to market on November 18, 2002
Olympus Guide Sheath	Olympus Medical Systems Corporation	K060243, cleared to market on June 23, 2006

5. Device Description

The LungPoint Tools are endoscopic tools used during bronchoscopy procedures. The LungPoint Sheath is designed to be used with a bronchoscope to provide a working

Broncus Medical, Inc.

Page 1 of 4

channel through which endoscopic tools, such as needles, dilation balloons, or other endoscopic devices may be introduced to the targeted lung tissue within the respiratory organs. The LungPoint Dilation Balloon is used to dilate tissue of the bronchial tree and may be inserted through the sheath or directly through the working channel of the bronchoscope. The bronchoscope is advanced to a predefined target following guidance of the LungPoint Software.

The materials used in the LungPoint Tools are commonly used medical grade materials and include platinum iridium markers.

Catalog Number	Working Length	Maximum Catheter Outer Diameter (OD)	Catheter Internal Diameter (ID)	Minimum accessory length for use through sheath
10007-1	900mm	2.65mm	2.0mm	965mm

The LungPoint Sheath has the following specifications:

The LungPoint Dilation Balloon has the following specifications:

Catalog Number	Balloon Size (OD x length)	Rated Balloon Pressure	Maximum Catheter Outer Diameter (OD)	Catheter Length/working length
10008-1	4mm x 6mm	20atm	1mm	1430mm/975mm

6. Intended Use

The LungPoint Tools are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs in a minimally invasive manner.

7. Comparison to Predicate Device

The LungPoint Tools are commonly used endoscopic tools with the same technological characteristics as the predicate devices. The indications for use of the LungPoint Tools for use with a bronchoscope guided by the LungPoint Software are all within the intended use of the predicate devices, which is to aid in reaching a target in the respiratory organ either directly as is the case with the sheath or through dilating target tissue with the dilation balloon. The technological characteristics of the subject devices are the same as those of the predicate devices with the following exceptions

• LungPoint Balloon: balloon size and length, balloon burst pressure and balloon catheter length as outlined.

	CRE Pulmonary Balloon Dilatation Catheter (K023337)	LungPoint Dilation Balloon
Balloon Size and Inflation	8mm @ 3 ATM	4mm @ 10 ATM
(balloon pressure)	9mm @ 5.5 ATM	
	10mm @ 9 ATM	
Catheter Length (cm)	155	143
Balloon Length (cm)	3.0	0.6
Rated Burst Pressure (atmospheres)	9	20

• LungPoint Sheath: stylet, catheter ID/OD and catheter length as outlined.

	Olympus Guide Sheath (K060243)	LungPoint Sheath
Stylet Provided	No	Yes - used to enhance pushability and to prevent airway mucosa entering sheath
Catheter ID/OD	2.1/2.7 mm	2.0/2.6 mm
Catheter Length	900 mm	975 mm

None of these differences raise new questions of safety and effectiveness. Performance of the subject devices has been verified by use of accepted methods.

8. Performance Data

The design and safety of the LungPoint Tools were verified by performing functional and performance testing. All tests were designed to subject the sheath and dilation balloon to stresses that exceed those which would be encountered during clinical use. Testing included the following:

- Dimensional testing
- Joint/tensile test
- Simulated use
- Balloon fatigue/burst pressure
- Balloon deflation time
- Radiopacity.

All testing results met the pre-determined acceptance criteria that were established in the test protocols. Based on the testing the LungPoint Tools are as safe, as effective, and perform at least as safely and effectively as the predicate devices.


9. Safety and Effectiveness

The LungPoint Tools labeling contains instructions for use and any necessary cautions and warnings, to provide for safe and effective use of the devices. The biocompatibility assessment of all patient contacting materials was performed in accordance with ISO 10993, *Biological Evaluation of Medical Devices*. Specifically, cytotoxicity, sensitization, intracutaneous reactivity and systemic toxicity (acute) were tested. In addition, the devices are sterilized using e-beam sterilization.

10. Conclusion

Based on the testing the LungPoint Tools are as safe, as effective, and perform at least as safely and effectively as the predicate devices.

Attachment 3: Predicate Device Labeling - Olympus Guide Sheath



JUN 2 3 2006

January 23 , 2006

SMDA 510(k) SUMMARY "Olympus Guide Sheath, XBO1-836-13"

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21CFR, Part 807, Subpart E, Section 807.92.

A. GENERAL INFORMATION

1. Applicant :	Olympus Medical Systems Corporation 2951 Ishikawa-cho, Hachioji-shi, Tokyo, 192-8507, Japan Registration number :8010047
2. Official Correspondent :	Laura Storm-Tyler Executive Director, Regulatory Affairs and Quality Assurance Olympus America Inc. Two Corporate Center Drive, Melville, NY 11747-9058, USA TEL 631-844-5688 FAX 631-844-5554 Registration Number :2429304
3. Manufacturer :	Aomori Olympus Co., Ltd. 2-248-1 Okkonoki Kuroishi-shi, Aomori-ken, 036-0357, Japan Registration Number : 9614641
B. DEVICE IDENTIFICATION	
1. Common/Usual Name :	Bronchoscope accessory
2. Device Name :	Olympus Guide Sheath, XBO1-836-13
3. Classification Name :	874.4680, class II,EOQ

C. PREDICATE DEVICES

Device Name	510(k) #	Manufacturer	Class	Product Code
Cytology Brush: BC-14/15/16C (EVIS 200 System)	#K931154	Olympus Corporation.	Π	EOQ

1 0 2

D. SUMMARY DESCRIPTION OF THE DEVICE

1. Summary

Olympus Medical Systems Corp., intends to introduce the Guide Sheath for use in the respiratory organs. This Guide Sheath,XBO1-836-13 has been designed to be used with Olympus bronchoscopes with 2.8 mm instrument channel, accessories and ultrasonic probe unit for performing diagnostic and therapeutic procedures. The Guide Sheath allows physicians to advance the endo-therapy accessories or guide the ultrasonic probe precisely to the targeted lesions repeatedly.

2. Materials

The patient contacting materials used as components of the Guide Sheath are identical to legally Marketed Olympus products.

E. SUMMARY

In summary, the Guide Sheath, XBO1-836-13 is basically identical to the predicate device in performance, materials and specifications. The subject device differs from the predicate device relative to the intended use.

F. INTENDED USE OF THE DEVICE

This instrument has been designed to be used with Olympus bronchoscopes, endo-therapy accessories and ultrasound probe to guide the endo-therapy accessories or ultrasound probe to the targeted area within the respiratory organs.

G. TECHNOLOGICAL CHARACTERISTICS

This Guide Sheath, XBO1-836-13 consists of tube sheath and piece of flouro tip, which shows the distal end position when X-ray was monitored. The following is the specification of the Guide Sheath;

Specifications	Dimension (m)
Outer Diameter	2.7
Inner Diameter	2.1
Sheath Length	900
Channel Size (Minimum ϕ)	2.8
Fluoro Tip for X-ray monitor	Provided

H. CLINICAL DATA

Olympus sponsored prospective clinical evaluation of the subject device which support the safe and effective use of the subject device for its intended use.

84



Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

JUN 2 3 2006

Olympus Medical Systems Corporation c/o Ms. Laura Storm-Tyler Two Corporate Center Dr. Melville, NY 11747-9058

Re: K060243

Trade/Device Name: Olympus Guide Sheath, XBO1-836-13 Regulation Number: 21 CFR 874.4680 Regulation Name: Bronchoscope and accessories (flexible or rigid) Regulation Class: II Product Code: EOQ Dated: May 30, 2006 Received: May 31, 2006

Dear Ms. Storm-Tyler:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Laura Storm-Tyler

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <u>http://www.fda.gov/cdrh/industry/support/index.html</u>.

Sincerely yours,

MB Esclemis, MD

Malvina B. Eydelman, M.D. Director Division of Ophthalmic and Ear, Nose and Throat Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

Indications for Use Statement

510(k) Number (if known) : K060243

Device Name: Olympus Guide Sheath, XBO1-836-13

Indications for Use :

This instrument has been designed to be used with Olympus bronchoscopes, endo-therapy accessories and ultrasound probe to guide the endo-therapy accessories or ultrasound probe to the target area within the respiratory organs.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use _____ (Per 21 CFR 801.109) OR

Over-The-Counter Use_____

(Optional Format 1-2-96)

(Division Sign-Off) Division of Ophthalmic Ear, Nose and Throat Devises

K060213 510(k) Number

7

Attachment 4: Instructions for Use LungPoint Dilation Balloon

Broncus Medical Inc.

LungPoint[™] Dilation Balloon

REF

Catalog number 10008


Instructions For Use (IFU)

Table of Contents

- 1 Device Description
- 2 Intended Use
- 3 Contraindications
- 4 Warnings / Precautions
- 5 Possible Complications
- 6 Preparation for Use
- 7 Operating Instructions
- 8 Storage

PN xxxxx_01 DRAFT





1 Device Description

The LungPoint Dilation Balloon is an endoscopic tool that is used to dilate the target lung tissue of the bronchial tree. The dilation balloon is inserted through the bronchoscope or other endoscopic tool, such as a sheath, and is used to dilate lung tissue. The balloon may be used during interventional bronchoscopy procedures as determined by the physician performing the bronchoscopy.

The LungPoint Dilation Balloon is comprised of a catheter with a non-compliant 4mm OD x 6mm long inflatable balloon at the distal end and a y-connector with female luer fittings at the proximal end. The side port of the Y connector is used for inflation and deflation of the balloon with a standard indeflator capable of 12atm pressure (typically with a gage) and a male luer lock connector. The center port contains a nitinol wire which is fixed in place for pushability and kink resistance during tortuous bronchoscopic navigation. Included are two radiopaque marker bands and the distal and proximal ends of the balloon inflation length for visualization of the device position under fluoroscopy

Product Configuration and Dimensions

The LungPoint Dilation Balloon (see illustration in Figure 1) has the following specifications:

Catalog	Balloon Size	Rated Balloon	Maximum Catheter	Catheter
Number	(OD x length)	Pressure	Outer Diameter (OD)	Length/working length
10008-1	4mm x 6mm	<u>20atm</u>	1mm	1430mm/975mm

Additionally, the LungPoint Dilation Balloon should be used with a flexible bronchoscope with a working channel of 2.0 mm or greater.

😽 Broncus Medical Inc.

2 Intended Use

The LungPoint Dilation Balloon is part of the LungPoint Tools, which are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs in a minimally invasive manner.

3 Contraindications

There are no known contraindications.




Figure 1: Dilation Balloon (shown inflated)

(1) WARNINGS AND PRECAUTIONS

4 Warnings / Precautions

- The LungPoint Dilation Balloon should only be used by physicians thoroughly trained in bronchoscopy and in the particular technique and procedure to be performed, and familiar with the associated risks. These operating instructions must be read, understood, and followed.
- The LungPoint Dilation Balloon should only be used in a manner consistent with the Bronchoscope manufacturer's Instructions for Use. Use of the LungPoint Dilation Balloon in contradiction to the instructions may harm the device, bronchoscope or patient.
- The rated burst pressure is 20 atmosphere, exceeding the rated burst pressure may result in harm to the patient.

Broncus Medical Inc.

5 Possible Complications

- Bleeding
- Infection

6 Preparation for Use

- 6.1 Materials Required
 - LungPoint Dilation Balloon
 - Indeflator capable of 12atm
 - Flexible bronchoscope with a working channel of 2.0 mm or greater
- 6.2 Device Inspection and Preparation
 - 6.2.1 Read all IFUs, package inserts, labels, and warnings for the LungPoint Dilation Balloon before beginning the procedure.

CAUTION: Inspect the sealed sterile device package before opening. If the seal is broken, contents may not be sterile and could pose a risk of patient infection.

6.2.2 Aseptically remove the LungPoint Dilation Balloon from the package and inspect for any damage, such as broken or crushed areas of the catheter shaft, sharp or protruding edges at the distal tip, or any kinks in the catheter shaft.

CAUTION: If any damage or irregularity is found, **DO NOT USE** the device as this could result in harm to the patient.

CAUTION: DO NOT USE the device past its "Use Before" (expiration) date.

CAUTION: DO NOT place or manipulate the catheter in the airways unless under direct visualization with a bronchoscope or using fluoroscopy as this may result in harm to the patient.

Broncus Medical Inc.

7 Operating Instructions

7.1 Carefully advance the balloon through the working channel of the bronchoscope using short (approximately 20 mm) increments, until the distal end of the catheter is within bronchoscopic view.

CAUTION: DO NOT force the catheter if resistance to insertion is encountered. Reduce the angulation of the bronchoscope until the device passes smoothly. This can result in (a) kinks to the device, (b) damage to the bronchoscope, or (c) harm to the patient such as punctures, hemorrhage or mucous membrane damage.

CAUTION: If any damage occurs, including a kink, DO NOT USE and discard.

7.2 Advance the catheter to the desired protrusion from the bronchoscope as required for navigation to the region of interest.

CAUTION: Do not angulate the bronchoscope abruptly while the device is extended from the distal end of the bronchoscope. This can result in kinks to the device or harm to the patient such as punctures, hemorrhage or mucous membrane damage.

- 7.3 Connect the indeflator via the luer lock connector and inflate to 10atm. Do not exceed the rated burst pressure of 20 atm.
- 7.4 If using fluoroscopy visualize balloon inflation location via the radiopaque marker bands and the distal and proximal ends of the balloon inflation length.
- 7.5 Deflate the balloon, and remove the device

CAUTION: If excessive resistance makes withdrawal difficult, reduce the angulation of the bronchoscope until the accessory can be withdrawn smoothly. Forcible withdrawal could damage the device or bronchoscope.

7.6 At the completion of the patient procedure, dispose of the device in accordance with applicable hospital, local, state and federal laws and regulations.

8 Storage

Store at controlled room temperature.

Broncus Medical Inc.

LEGAL NOTICE

© 2013 Broncus Medical, Inc. All rights reserved.

Broncus, the Broncus logo, and LungPoint are trademarks or registered trademarks of Broncus Medical, Inc.

Patents Pending.

REF	Catalog number
Ronly	Federal law (USA) restricts this device to sale by or on the order of a physician.
<u></u>	Caution: consult instructions for use
STERILE R	Irradiation sterilized. Sterility guaranteed if package unopened and undamaged.
LOT	Serial number or batch code
	Use before
2	For single use
	Manufacturer

Graphic Symbol Legend for Medical Device Labeling



Broncus Medical Inc.

Manufacturer: Broncus Medical, Inc. 1400 N. Shoreline Blvd., Bldg. A, Suite 8 Mountain View, CA, 94043 USA (877) 428-1600

Attachment 5: Instructions for Use LungPoint Sheath

Broncus Medical Inc.

LungPoint[™] Sheath

REF


Catalog number 10007

Instructions For Use (IFU)

Table of Contents

- 1 Device Description
- 2 Intended Use
- 3 Contraindications
- 4 Warnings / Precautions
- 5 **Possible Complications**
- 6 Preparation for Use
- 7 Operating Instructions
- 8 Storage

PN xxxxx_01 DRAFT



1 Device Description

The LungPoint Sheath is an endoscopic tool that is designed to be used with bronchoscopes to provide a working channel through which endoscopic tools, such as needles, dilation balloons, or other endoscopic devices may be introduced to the targeted area within the respiratory organs.

The sheath enables physicians to easily access targets and allows for multiple approaches to the preselected target. It allows for the repeated placement of endoscopic tools to a specified lesion(s) during one procedure. The sheath is provided with a stylet, which is used to minimize the amount of airway mucosa entering the sheath's lumen and/or to provide rigidity (i.e. pushability).

The sheath (with stylet, if necessary) is inserted through a standard bronchoscope. It is then advanced to a target under the guidance of the LungPoint Software. Placement of the sheath may require the use of other endoscopic tools. Once placed, other endoscopic tools, such as needles, dilation balloons, or other endoscopic devices, may be introduced through the sheath to the targeted area within the respiratory organs. If a stylet is used, it must be removed prior to the introduction of the other tools.

The LungPoint Sheath is comprised of a braid reinforced tubing, to resist kinking with articulation along with an incorporated stylet. The stylet when mated, provides for a rounded tip during insertion and pushability for navigation to the target. Removal of the stylet allows for standard 2.0mm working channel bronchoscopic accessories to be used though the lumen of the sheath. The tip of the Sheath is marked visually with black bands at the tip and in 10mm increments to 60mm to provide the user with an indication of the traveled depth during insertion. Also included are multiple radiopaque marker bands at the distal end at the tip, 5mm, 10mm and 20mm to aid visualization of the sheath under fluoroscopy.

Broncus Medical Inc.

Product Configuration and Dimensions

The LungPoint Sheath (see illustration in Figure 1) has the following specifications:

Catalog Number	Working Length	Maximum Catheter Outer Diameter (OD)	Catheter Internal Diameter (ID)	Minimum accessory length for use through sheath
10007-1	900mm	2.65mm	2.0mm	965mm

Additionally, the LungPoint Sheath should be used with a flexible bronchoscope with a working channel of 2.8 mm or greater.




Figure 1: LungPoint Sheath

2 Intended Use

The LungPoint Sheath is part of the LungPoint Tools, which are endoscopic tools used with bronchoscopes and intended to be used as accessories to the LungPoint Software to aid in reaching a targeted area within the respiratory organs in a minimally invasive manner.

3 Contraindications

There are no known contraindications.

В

Broncus Medical Inc.

WARNINGS AND PRECAUTIONS

4 Warnings / Precautions

- The LungPoint Sheath should only be used by physicians thoroughly trained in bronchoscopy and in the particular technique and procedure to be performed, and familiar with the associated risks. These operating instructions must be read, understood, and followed.
- DO NOT activate active devices (e.g., electrosurgical devices) inside the lumen of the sheath. The active portion of the device must be extended beyond the sheath prior to activation.
- While advancing the sheath, DO NOT USE excessive force as this may advance the device in an uncontrolled manner that could result in harm to the patient's airway wall and cause bleeding.
- The LungPoint Sheath should only be used in a manner consistent with the Bronchoscope manufacturer's Instructions for Use. Use of the LungPoint Sheath in contradiction to the instructions may harm the device, bronchoscope or patient.

5 Possible Complications

- Bleeding
- Infection

6 Preparation for Use

6.1 Materials Required

- LungPoint Sheath
- Luer-lock aspiration syringe
- Flexible bronchoscope with a working channel of 2.8 mm or greater

6.2 Device Inspection and Preparation

6.2.1 Read all IFUs, package inserts, labels, and warnings for the device before beginning the procedure.

CAUTION: Inspect the sealed sterile device package before opening. If the seal is broken, contents may not be sterile and could pose a risk of patient infection.

6.2.2 Aseptically remove the LungPoint Sheath from the package and inspect for any damage, such as broken or crushed areas of the catheter shaft, sharp or protruding edges at the distal tip, or any kinks in the catheter shaft.

Broncus Medical Inc.

CAUTION: If any damage or irregularity is found, DO NOT USE the device as this could result in harm to the patient.

CAUTION: DO NOT USE the device past its "Use Before" (expiration) date.

6.2.3 Before inserting the device into the working channel of the bronchoscope, verify that the stylet is secured to the sheath via the luer lock connection.

7 Operating Instructions

7.1 Carefully advance the Sheath through the working channel of the bronchoscope using short (approximately 20 mm) increments, until the distal end of the catheter is within bronchoscopic view.

CAUTION: DO NOT force the catheter if resistance to insertion is encountered. Reduce the angulation of the bronchoscope until the device passes smoothly. This can result in (a) kinks to the device, (b) damage to the bronchoscope, or (c) harm to the patient such as punctures, hemorrhage or mucous membrane damage.

CAUTION: If any damage occurs, including a kink, DO NOT USE and discard.

7.2 Advance the catheter to the desired protrusion from the bronchoscope as required for navigation to the region of interest.

CAUTION: Do not angulate the bronchoscope abruptly while the device is extended from the distal end of the bronchoscope. This can result in kinks to the device or harm to the patient such as punctures, hemorrhage or mucous membrane damage.

7.3 For accessory use: remove the stylet by loosening the luer lock connection and withdraw the stylet. Insert the accessory device into the lumen.

CAUTION: If excessive resistance makes withdrawal difficult, reduce the angulation of the bronchoscope until the accessory can be withdrawn smoothly. Forcible withdrawal could damage the sheath, accessory or bronchoscope.

- 7.4 If using fluoroscopy visualize radiopaque accessory devices position relative to the tip marker band.
- 7.5 Withdraw the Sheath from the bronchoscope.
- 7.6 At the completion of the patient procedure, dispose of the device in accordance with applicable hospital, local, state and federal laws and regulations.

8 Storage

Store at controlled room temperature.

Broncus Medical Inc.

© 2013 Broncus Medical, Inc. All rights reserved.

Broncus, the Broncus logo, and LungPoint are trademarks or registered trademarks of Broncus Medical, Inc.

Patents Pending.

LEGAL NOTICE

REF	Catalog number
Ronly	Federal law (USA) restricts this device to sale by or on the order of a physician.
	Caution: consult instructions for use
STERILE R	Irradiation sterilized. Sterility guaranteed if package unopened and undamaged.
LOT	Serial number or batch code
	Use before
2	For single use
	Manufacturer

Graphic Symbol Legend for Medical Device Labeling



Broncus Medical Inc.

Manufacturer: Broncus Medical, Inc. 1400 N. Shoreline Blvd., Bldg. A, Suite 8 Mountain View, CA, 94043 USA (877) 428-1600 **Attachment 6: Literature Reference**

Dimensions of the Normal Human Trachea

Éamann Breatnach¹ Gypsy C. Abbott² Robert G. Fraser¹ The coronal and sagittal diameters of the tracheal air column were measured on posteroanterior and lateral chest radiographs of 808 patients with no clinical or radiographic evidence of respiratory disease. The 430 male and 378 female subjects were 10–79 years of age. Assuming a normative range that encompasses three standard deviations from the mean or 99.7% of the normal population, the upper limits of normal for coronal and sagittal diameters, respectively, in men aged 20–79, are 25 mm and 27 mm; in women, they are 21 mm and 23 mm, respectively. The lower limit of normal for both dimensions is 13 mm in men and 10 mm in women. Deviation from these figures reflects pathologic widening or narrowing of the tracheal air column. No statistically significant correlation was found between tracheal caliber and body weight or body height.

Changes in tracheal dimensions occur in a variety of conditions. For example, generalized widening is a characteristic feature of tracheobronchomegaly and tracheomalacia; generalized narrowing is seen in tracheobronchopathia osteochondroplastica and may be a feature of relapsing polychondritis. Knowledge of normal tracheal dimensions on conventional chest radiographs is essential to the diagnosis of these conditions. Previously published figures for normal tracheal caliber show considerable variation and generally are based on sample populations of insufficient size to establish statistical validity for all decades. We measured the coronal and sagittal tracheal diameters on posteroanterior and lateral chest radiographs of a large sample of patients with no known respiratory disease to establish normal ranges for these measurements in men and women of various ages and to determine whether correlation exists between tracheal width and body weight or body height.

Subjects and Methods

The chest radiographs of 808 patients (430 males and 378 females) were examined. Radiographs were obtained on maximal inspiration in posteroanterior and lateral projections using a focus-film distance of 10 feet (3.05 m) and a 6-inch (15.2 cm) air gap for scatter reduction. This technique results in a mid-plane magnification factor of 1.08, slightly less than the magnification factor of 1.10 obtained with the more conventional 6-foot (1.83 m) focus-film distance and grid technique. Exposures were made at 145 kVp and 500 mA for the posteroanterior projection and at 110 kVp and 700 mA for the lateral projection. All patients were ambulatory and had neither historical nor radiographic evidence of respiratory disease. Patients were included in the study only when their chest radiographs on maximal inspiration were technically adequate for accurate measurement of the intrathoracic trachea in both projections, and in whom local or general tracheal deformity from any cause was absent. Patients' age, height, and weight were recorded under supervision at the time of the radiographic examination. The internal diameter of the tracheal air column was measured at a level 2 cm above the projected top of the aortic arch on both posteroanterior and lateral radiographs to obtain coronal and sagittal dimensions, respectively (fig. 1). Measurements

Received July 6, 1983; accepted after revision December 14, 1983.

¹ Department of Diagnostic Radiology, University of Alabama, Birmingham Medical Center, 619 S. 19th St., Birmingham, AL 35233. Address reprint requests to E. Breatnach.

² Department of Statistics, University of Alabama, Birmingham Medical Center, Birmingham, AL 35233.

AJR 141:903–906, May 1984 0361–803X/84/1425–0903 © American Roentgen Ray Society 904

BREATNACH ET AL.




Fig. 1.—Level of measurement (solid line) of coronal (A) and sagittal (B) tracheal diameters, 2 cm above projected top of aortic arch (broken line).

were made of the air column alone, excluding the tracheal wall. This level was chosen because of the high visibility of the aortic arch on diagnostic chest radiographs, whereas the carina is not always precisely identifiable. Internal tracheal diameters were used for all measurements.

Results

Tracheal dimensions were grouped according to the subjects' age and gender, and the mean coronal and sagittal diameters (±SD) were calculated for each decade for males and females (table 1). In men, at all ages from 20 to 79 years, statistically significant differences ($\rho = 0.05$) were observed between sagittal and coronal diameters, the former being consistently greater than the latter (figs. 2 and 3). In women, although mean sagittal diameters were greater than coronal from the third decade on, the difference was not as great as in men and assumed statistical significance (p = 0.05) only in subjects over 49 years of age. In the third through eighth decades of life, the coronal and sagittal dimensions were greater in men than in women; the differences were statistically significant (p = 0.05). During the second decade, although the mean values in both dimensions were greater in boys than in girls, the differences between the genders were not statistically significant (fig. 4). Gender differences were observed in growth patterns (figs. 2 and 3): In men, both dimensions increased rather sharply during the second and third decades, then more slowly (but significantly) during the fourth decade; in women, growth was much more gradual from the second through fourth decades. In neither gender did any change occur in either dimension from the fourth decade on.

Relations between tracheal measurements and body weight and height were analyzed by multiple regression techniques, including cross-validation; no statistically significant correlation was found.

Discussion

The most commonly quoted figures for tracheal caliber in normal subjects are those of Katz et al. [1], Jesseph and Merendino [2], Greene [3], and Fraser and Paré [4]; data from these four studies are summarized in table 2. Although values for normal ranges are quoted elsewhere [5, 6], the data on which the figures are based generally lack statistical validity. It can be seen from table 2 that there is considerable variation in the normal ranges published. Greene [3] measured the internal tracheal diameter 1 cm above the aortic arch, whereas Jesseph and Merendino [2] used a point just proximal to the tracheal bifurcation. We found a point 2 cm above the aortic arch to be more consistently reproducible, though we observed no real variation between measurements taken at either level. The reported differences in normal ranges probably reflect, at least in part, differences in age groups and gender predominance among the populations studied. In the present study, we used a larger group of patients with convincing numeric representation for both genders in each decade from 10 to 79 years of age. The small size of the SD for each mean in our series indicates satisfactory numeric representation for each decade on a statistical basis and the validity of comparative data between the genders. The increase in coronal and sagittal tracheal dimensions for both genders up to age 30 (fig. 4) is important in assessing normal values for younger patients. The similarity between coronal and sagittal diameters in the second decade in males and in the second and third decades in females correlates well with the work of Engle [7], who found that the trachea in children and young adults is roughly circular in cross section. In all decades after ages 20 (men) and 30 (women), the trachea assumes a more ovoid shape in cross section; the differences betwen coronal and sagittal diameters are greater in men than in women.

Changes in tracheal dimensions occur in a variety of conditions; for example, there is dilatation in tracheobronchomegaly [1, 8] (fig. 5) and tracheomalacia [9] and narrowing in tracheobronchopathia osteochondroplastica [10], relapsing polychondritis [11], fibrosing mediastinitis [12, 13], rhinoscleroma [14], and "saber sheath" trachea [15]. The lastnamed deserves brief description.

Saber-sheath trachea, as defined by Greene and Lechner [15], involves a coronal tracheal diameter less than two-thirds the sagittal diameter (fig. 6). In Greene's later report [3], the mean coronal and sagittal diameters in a control group (60 men; mean age, 66.4 years) were 19.7 and 22.5 mm, respectively, and in a study group (60 men; mean age, 63.9 years) 12.0 and 25.8 mm, respectively. Such a pronounced difference between the two diameters was not seen in any of our patients, indicating that saber-sheath trachea is not a natural consequence of aging, as was suggested by Simmonds [16], cited by Greene [3]. That the deformity does not occur in the normal population at any age adds support to Greene's

AJR:142, May 1984

NORMAL TRACHEAL DIMENSIONS

		Male			Female	
Age Group (years)		Tracheal Diameter	(mm): Mean ± SD	No. of Cubicato	Tracheal Diameter	(mm): Mean ± SD
	No. of Subjects	Coronal	Sagittal	No. of Subjects	Coronal	Sagittal
10–19	26	15.5 ± 2.8	15.4 ± 3:1	22	14.4 ± 1.6	14.5 ± 1.3
20–29	81	18.7 ± 2.0	19.3 ± 2.0	98	15.7 ± 1.6	15.6 ± 1.7
30–39	72	19.2 ± 2.1	19.7 ± 2.4	64	16.0 ± 1.8	16.3 ± 2.3
40–49	69	19.5 ± 2.3	20.3 ± 2.2	45	16.6 ± 2.0	16.8 ± 2.2
50–59	80	19.2 ± 2.3	20.4 ± 2.6	65	16.5 ± 1.6	17.0 ± 2.0
60–69	71	19.5 ± 2.2	20.7 ± 2.5	48	16.8 ± 2.0	17.2 ± 2.3
70–79	31	19.7 ± 2.2	20.8 ± 1.8	36	16.4 ± 2.4	16.5 ± 2.3

Note.—Coronal and sagittal tracheal diameters were defined as the internal diameters of the tracheal air column as measured at a level 2 cm above the projected top of the aortic arch on posteroanterior and lateral radiographs, respectively.




Fig. 2.—Mean sagittal tracheal diameters (mm) in normal subjects by gender and age group (decades).




Fig. 3.—Mean coronal tracheal diameters (mm) in normal subjects by gender and age group (decades).

contention [15] that it is almost invariably associated with chronic obstructive airway disease: Of the 60 patients with saber-sheath trachea, 95% had clinical evidence of chronic obstructive pulmonary disease, as compared with only 18% in the control group. The value of this radiographic sign is twofold: It indicates the presence of chronic obstructive pul-




Fig. 4.—Mean coronal tracheal diameters (mm) in normal subjects by gender and age (years) during second decade of life.

monary disease when other convincing radiographic evidence is lacking, and it militates against the erroneous assumption that coronal narrowing is caused by a mediastinal mass.

In our study, no correlation was observed between tracheal dimensions and subjects' weight or height. This contrasts with the results reported by Fraser and Paré [4] in a study of 350 chest radiographs of roughly equal numbers of asymptomatic men and women 30–80 years of age; these authors found that sagittal tracheal diameters increased gradually with increasing height, although neither age nor weight had any apparent influence. The average sagittal diameter increased from 17.6 mm in subjects 50–59 inches (127–150 cm) tall to 22.5 mm in subjects 70–72 inches (178–183 cm) tall; corresponding coronal dimensions were not recorded. We have no explanation for the discrepancy between the two studies regarding the effect of subjects' height on tracheal dimensions.

Assuming a normative range that encompasses three SDs from the mean, the upper limits of normal for coronal and sagittal diameters of the tracheal air column in men aged 20–79 are 25 and 27 mm, respectively; in women aged 20–79, they are 21 and 23 mm, respectively. The lower limit of normal for both dimensions is 13 mm in men and 10 mm in women. Any deviation from these figures indicates pathologic widening or narrowing of the tracheal air column.

ACKNOWLEDGMENT

We thank Sharon Compton for assistance with computer analysis.

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 30 Attech ment 6 Page 4 of 5

BREATNACH ET AL.

AJR:142, May 1984

TABLE 2: Normal Coronal and Sagittal Tracheal Diameters Reported by Previous Investigators

Reference	No. of Subjects	Age Range (mean)	Tracheal Diameter (mm): Range (mean ± SD)		
	(gender)	in Years	Coronal	Sagittal	
[1]	50		13-25 (20.2 ± 3.4)		
[2]	21 (M) 26 (F)	13-86	15-27 (22.0 ± 2.6) 13-25 (17.0 ± 2.2)		
[3]	60 (M)	(66.4)	15-26 (19.7 ± 2.0)	18-32 (22.5 ± 2.4)	
[4]	175 (M) 175 (F)	30-80	11–26 (19.5) 9–24 (15.5)	11–30 (17.5) 7–26 (15.5)	




Fig. 5.—Tracheobronchomegaly (Mounier-Kuhn syndrome). Coronal diameter (A, solid line) = 31 mm; sagittal diameter (B, solid line) = 22 mm. (Broken line marks level of top of aortic arch.)

REFERENCES

- 1. Katz I, LeVine M, Herman P. Tracheobronchiomegaly: the Mounier-Kuhn syndrome. AJR 1962;88 : 1084–1094
- 2. Jesseph JE, Merendino KA. Dimensional interrelationships of the major components of the human tracheobronchial tree. Surg Gynecol Obstet 1957;105 : 210-214
- 3. Greene R. "Saber-sheath" trachea: relation to chronic obstructive pulmonary disease. AJR 1978;130 : 441-445
- 4. Fraser RG, Paré JAP. Diagnosis of diseases of the chest, vol 1. Philadelphia: Saunders, 1977 : 56
- Schinz HR, Baensch WE, Frommhold W, Glauner R, Uehlinger 5. E, Wellauer J. Roentgen diagnosis, 2d ed, vol 4, pt 2. New York: Grune & Stratton, 1973 : 1
- 6. Shanks SC, Kerley P. A textbook of x-ray diagnosis by British authors, 4th ed, vol 3. Philadelphia: Saunders, 1973 : 88
- 7. Engle S. The child's lung. Ann Arbor: University Microfilms, 1955
- 8. Mounier-Kuhn P. Dilatation de la trachée; constatations radiographiques et bronchoscopiques. Lyon Med 1932;150 : 106-109




Fig. 6.--Saber-sheath trachea. Coronal diameter (A, solid line) = 16 mm; sagittal diameter (B, solid line) = 30 mm. (Broken line marks level of top of aortic arch.)

- 9. Johnson TH, Mikita JJ, Wilson RJ, Feist JH. Acquired tracheomalacia. Radiology 1973;109 : 577-580
- 10. Howland WJ Jr, Good CA. The radiographic features of tracheopathia osteoblastica. Radiology 1958;71: 847-850
- 11. Gobson GJ, Davis P. Respiratory complications of relapsing polychondritis. Thorax 1974;29: 726-731
- 12. Schowengerdt CG, Suyemoto R, Main FB. Granulomatous and fibrous mediastinitis. A review and analysis of 180 cases. J Thoracic Cardiovasc Surg 1969;57: 365-379
- 13. Barrett NR. Idiopathic mediastinal fibrosis. Br J Surg 1958; 46:207-218
- 14. Miller RH, Shulman JB, Canalis RF, et al. Klebsiella rhinoscleromatis: a clinical and pathogenic enigma. Otolaryngol Head Neck Surg 1979;87 : 212-221
- 15. Greene R, Lechner GL. "Saber-sheath" trachea: a clinical and functional study of marked coronal narrowing of the intrathoracic trachea. Radiology 1975;115 : 265, 268
- 16. Simmonds M. Uber Alterssabelscheidentrachea. Virchows Arch [Pathol Anat] 1905;179 : 15-28

Attachment 7: LungPoint Balloon, Simulated Use Testing including Deflation Time Test Data

Attachment 8: Dilation Balloon Drawings (Revision A and B)

Attachment 9: LungPoint Tools Shelf-life Report

Records processed under FOIA Request # 2015-7253; Released by CDRH on 02-01-2016 34/5/

FDA CDRH DMC

AUG 27 2013

Date: August 26, 2013

Received

U.S. Food and Drug Administration Center for Devices and Radiological Heath Document Mail Center – WO66-0609 10903 New Hampshire Avenue Silver Spring, MD 20993-0002

RE: Premarket Notification 510(k) Electronic Copy:

K131234 Broncus Medical, Inc. LungPoint Tools

To Whom It May Concern:

The electronic files (eCopy) are an exact duplicate of the paper copy included with this submission.

If you should have any questions regarding this submission please contact me at (b)(4) or email or email (b)(4) or email (b)(4) or email (b)(4) (b)(4)

Sincerely,

Responsible Third Party Official



D)(4) - 30 May 03 Page 1 of 1

Revis

ID at CDRH-FOISTATUS@fda.hhs.gov or \$01-796-8118

Date: August 26, 2013

U.S. Food and Drug Administration Center for Devices and Radiological Heath Document Mail Center – WO66-0609 10903 New Hampshire Avenue Silver Spring, MD 20993-0002

RE: Additional Information for

K131234 Broncus Medical, Inc. LungPoint Tools

To Whom It May Concern:

Enclosed in duplicate is the additional information received from Broncus Medical regarding K131234. This information has been reviewed in detail (see attached) review addendum and was found to answer the questions raised from a letter dated June 12, 2013 from Srinivas "Nandu" Nandkumar, Ph.D..

Recommendation: Substantial Equivalence

If you should	d have any fu	rther questions	regarding th	is submission	please	contact
me at (b)(4)	or fax	(b)(4)	or email at	(b)(4)		Please
fax any cor	respondence	regarding this	submission	tc <mark>(b)(4)</mark>		
(b)(4)	ļ					

	Sincerely
(b)(6)	
	Responsible Third Party Official

(b)(4) Revision Effective 16 Feb 04 r 301 796 8118 Page 1 of 1

t CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118