

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
510(k) Summary

---

FEB 21 2014

## 5 510(k) Summary

### 5.1 Submitter's Contact Information

Straumann USA, LLC (on behalf of Institut Straumann AG)

60 Minuteman Road

Andover, MA 01810

Phone Number: 1-978-747-2509

Fax Number: 1-978-747-0023

Contact Person: Jennifer M. Jackson, MS

Date of Submission: 20-Feb-2014

### 5.2 Name of the Device

Trade Name: Straumann® Variobase™ Abutments

Common Name: Dental Implant Abutment

Classification Name: Abutment, Implant, Dental, Endosseous

Regulation Number: §872.3630

### 5.3 Predicate Device(s)

- K120822, Straumann CARES Variobase Abutment NNC, RN, WN, NC, RC (Institut Straumann AG)
- K111935, Ti-Base Abutment (NT-Trading GmbH & Co. KG)

### 5.4 Device Description

The Straumann® Variobase™ Abutments are pre-manufactured (stock) abutments, sometimes referred to as "Ti-bases". Straumann® Variobase™ Abutments are available to fit Straumann® dental implant platforms NNC (Narrow Neck CrossFit®), RN (Regular Neck), WN (Wide Neck), NC (Narrow CrossFit®), and RC (Regular CrossFit®). A dental laboratory technician would design the corresponding coping and/or crown (the second component of the Variobase two-piece abutment) and/or prosthetic restoration in the



# **Traditional 510(k) Submission**

## **Straumann® Variobase™ Abutments**

### 510(k) Summary

---

dental laboratory using either a burnout coping or STL model for open CAD software. The coping and/or crown would be manufactured via validated Straumann milling.

#### **5.5 Intended Use**

Straumann® Variobase™ Abutments are pre-manufactured (stock) abutments intended to be placed onto Straumann dental implants to provide support for customized prosthetic restorations such as crowns and bridges. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.

#### **5.6 Technological Characteristics**

Straumann® Variobase™ Abutments are two-piece abutments consisting of a pre-manufactured (stock) abutment made from a titanium-aluminum-niobium alloy and a coping and/or crown which is designed in the dental laboratory by a dental technician using open CAD software and manufactured via validated Straumann milling.

The Ti-base components of the Straumann® Variobase™ Abutments are identical to the Ti-base components of the Straumann predicate (K120822). The Ti-base components are also equivalent to the Ti-base components identified in K111935.

The materials which may be used to manufacture the coping/crown component of the Straumann® Variobase™ Abutments are identical to the identified predicate devices and include:

Milling:            Polycon® ae (temporary restorations – K120822)  
                          Zerion® (K120822)

#### **5.7 Performance Testing**

The material used in the manufacture of Straumann® Variobase™ Abutments is a titanium-aluminum-niobium alloy which meets the requirements of ISO 5832-11. Bench testing was performed with Polycon® ae and Zerion® to evaluate the fatigue load limits of the proposed Straumann® Variobase™ Abutments. Dynamic fatigue tests were conducted in accordance to the FDA guidance document *“Guidance for Industry and FDA*

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

510(k) Summary

---

*Staff – Class II Special Controls Guidance Document: Root-form Endosseous Dental Implants and Endosseous Dental Abutments”.*

**5.8 Conclusion**

The documentation submitted in this premarket notification demonstrates that the Straumann® Variobase™ Abutments are substantially equivalent to the predicate devices and do not pose new issues of safety and effectiveness when used as labeled.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

February 21, 2014

Straumann USA, Limited Liability Company  
Jennifer M. Jackson, MS  
60 Minuteman Road  
Andover, MA 01810

Re: K132219

Trade/Device Name: Straumann® Variobase™ Abutments  
Regulation Number: 21 CFR 872.3630  
Regulation Name: Abutment, Implant, Dental, Endosseous  
Regulatory Class: II  
Product Code: NHA  
Dated: January 23, 2014  
Received: January 24, 2014

Dear Ms. Jackson:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Kwame O.  
Ulmer-S**

for

Erin I. Keith, M.S.  
Acting Director  
Division of Anesthesiology, General Hospital,  
Respiratory, Infection Control and  
Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): **K132219**

Device Name: **Straumann® Variobase™ Abutments**

### Indications for Use:

The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.

Prescription Use   **X**    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

The FDA logo is centered at the bottom of the page. It features the letters 'FDA' in a stylized, bold font. Overlaid on the logo is handwritten text in black ink. The text includes the name 'Mary A. Runner-S', the device name 'Straumann', the 510(k) number 'K132219', and a date '12/22/05'.

Page 1 of 1



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

February 21, 2014

Straumann USA, Limited Liability Company  
Jennifer M. Jackson, MS  
60 Minuteman Road  
Andover, MA 01810

Re: K132219

Trade/Device Name: Straumann® Variobase™ Abutments  
Regulation Number: 21 CFR 872.3630  
Regulation Name: Abutment, Implant, Dental, Endosseous  
Regulatory Class: II  
Product Code: NHA  
Dated: January 23, 2014  
Received: January 24, 2014

Dear Ms. Jackson:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Kwame O.**  
**Ulmer-S**



for

Erin I. Keith, M.S.  
Acting Director  
Division of Anesthesiology, General Hospital,  
Respiratory, Infection Control and  
Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Concurrence & Template History Page**  
[THIS PAGE IS INCLUDED IN IMAGE COPY ONLY]

Full Submission Number: K132219

For Office of Compliance Contact Information:

[http://insideportlets.fda.gov:9010/portal/page?\\_pageid=197,415881&\\_dad=portal&\\_schema=PORTAL&org=318](http://insideportlets.fda.gov:9010/portal/page?_pageid=197,415881&_dad=portal&_schema=PORTAL&org=318)

For Office of Surveillance and Biometrics Contact Information:

[http://insideportlets.fda.gov:9010/portal/page?\\_pageid=197,415881&\\_dad=portal&\\_schema=PORTAL&org=423](http://insideportlets.fda.gov:9010/portal/page?_pageid=197,415881&_dad=portal&_schema=PORTAL&org=423)

Digital Signature Concurrence Table	
Reviewer Sign-Off	Michael Mendelson
Branch Chief Sign-Off	Susan Runner
Division Sign-Off	Kwame O. Ulmer - S 2014.02.21 14:00:59 -05'00'

Template Name: K1(A) – SE after 1996

(b)(4) Trade Secret Process





## Indications for Use

510(k) Number (if known): **K132219**

Device Name: **Straumann® Variobase™ Abutments**

### Indications for Use:

The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.

Prescription Use   **X**    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF  
NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)



Page 1 of 1



K132219

FDA CDRH DMC

July 16, 2013

JUL 17 2013

U.S. Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center – W066-G609  
10903 New Hampshire Avenue  
Silver Spring, Maryland 20993-0002

Received

**Subject: Traditional 510(k) Premarket Notification: Straumann® Variobase™ Abutments**

Dear Sir or Madam:

On behalf of Institut Straumann AG, Straumann USA, LLC submits this Traditional 510(k) Premarket Notification for the Straumann® Variobase™ Abutments in accordance with 21 CFR 807.81. The intended use and fundamental operating principles of the proposed devices are substantially equivalent to previously cleared devices as detailed in this premarket submission.

**Submitter:**

Straumann USA, LLC (on behalf of Institut Straumann AG)  
60 Minuteman Road  
Andover, MA 01810

**Primary Contact:**

Jennifer M. Jackson, MS  
Senior Regulatory Affairs Project Manager  
Straumann USA, LLC  
60 Minuteman Road  
Andover, MA 01810  
Telephone: 800-448-8168 x2509  
Fax: 978-747-0023

**Classification Name of Device:**

Classification Name:	Abutment, Implant, Dental, Endosseous
Device Product Code:	NHA
Product Classification:	Class II
Panel:	Dental
Regulation Number:	§872.3630

8



**Design and Use of the Device:**

Question	YES	NO
Is the device intended for prescription use (21 CFR 801 Subpart D)?	X	
Is the device intended for over-the counter use (21 CFR 807 Subpart C)?		X
Does the device contain components derived from a tissue or other biologic source?		X
Is the device provided sterile?		X
Is the device intended for single use?	X	
Is the device a reprocessed single use device?		X
If yes, does this device type require reprocessed validation data?		
Does the device contain a drug?		X
Does the device contain a biologic?		X
Does the device use software?		X
Does the submission include clinical information?		X
Is the device implanted?		X

The term “substantially equivalent” as used herein is intended to be a determination of substantial equivalence from an FDA-regulatory point of view under the Federal Food, Drug and Cosmetic Act and relates to the fact that the product can be marketed without premarket approval or reclassification. These products may be considered distinct from a patent point of view. The term “substantially equivalent” is not applicable to, and does not diminish, any patent claim related to this product or the technology used to manufacture the product.

Information in this premarket notification is considered proprietary or trade secret or confidential commercial information. The company requests that all such information not be disclosed pursuant to 18 U.S.C. §1905, 5 U.S.C. §552, 21 U.S.C. §331(j), and all other applicable laws and regulations.

Two copies of the Traditional 510(k) Premarket Notification are enclosed. The second copy is being provided in PDF format on CD and is an exact duplicate of the hardcopy. Further, in accordance with the Medical Device User Fee and Modernization Act of 2001 (“MDUFMA”), Straumann USA, LLC has submitted the appropriate application fees. A copy of the User Fee Cover Sheet is provided with the enclosed 510(k) Premarket Notification.

We trust that the foregoing information will be sufficient to permit FDA to make a finding of substantial equivalence for the proposed Straumann® Variobase™ Abutments to the currently marketed devices as presented in this premarket notification.

Please address any questions regarding this 510(k) Premarket Notification to the undersigned.

Sincerely,

Jennifer M. Jackson, MS  
Senior Regulatory Affairs Project Manager

Enclosures

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

---











### Table of Contents

1	Medical Device User Fee Cover Sheet (Form FDA 3601) .....	1-1
2	CDRH Premarket Review Submission Cover Sheet .....	2-1
3	510 (k) Cover Letter.....	3-1
4	Indications for Use Statement.....	4-1
5	510(k) Summary.....	5-1
5.1	Submitter’s Contact Information .....	5-1
5.2	Name of the Device.....	5-1
5.3	Predicate Device(s).....	5-1
5.4	Device Description .....	5-1
5.5	Intended Use.....	5-1
5.6	Technological Characteristics.....	5-2
5.7	Performance Testing.....	5-2
5.8	Conclusion .....	5-2
6	Truthful and Accuracy Statement.....	6-1
7	Class III Summary and Certification.....	7-1
8	Financial Certification or Disclosure Statement.....	8-1
9	Declarations of Conformity and Summary Reports .....	9-1
9.1	ISO 10993-1:2009 .....	9-3
9.2	ISO 10993-5:2009 .....	9-5
9.3	ISO 10993-12:2007 .....	9-7
9.4	ISO 10993-18:2005 .....	9-9
9.5	ISO 5832-11:1994 .....	9-11
9.6	ISO 17664:2004 .....	9-13
9.7	ISO 17665-1:2006 .....	9-15

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

---

9.8	ISO/TS 17665-2:2009 .....	9-17
9.9	ISO 14801:2007 .....	9-19
10	Device Description .....	10-1
		
		
		
		
		
		
		
11	Substantial Equivalence Discussion.....	11-1
12	Proposed Labeling .....	12-1
12.1	Package Label .....	12-1
12.2	Proposed Package Insert/Instructions for Use for Straumann® Variobase™ Abutments.....	12-2
13	Sterilization and Shelf Life.....	13-1
13.1	Sterilization.....	13-1
13.2	Shelf Life .....	13-2
13.3	Packaging .....	13-3
14	Biocompatibility .....	14-1
		
15	Software .....	15-1
16	Electromagnetic Compatibility and Electrical Safety.....	16-1
17	Performance Testing – Bench.....	17-1
		
		

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process	[Redacted]
ss	[Redacted]
[Redacted]	[Redacted]
[Redacted]	[Redacted]
[Redacted]	[Redacted]
[Redacted]	[Redacted]
[Redacted]	[Redacted]

18	Performance Testing – Animal .....	18-11
19	Performance Testing – Clinical .....	19-1
	Appendix 1 – Engineering Drawings .....	A1-1
	Appendix 2 – NT-Trading Product Catalog .....	A2-1
	Appendix 3 – NT-Trading Ti-Base Package Insert .....	A3-1
	Appendix 4 – NT-Trading Ti-Base 510(k) Summary .....	A4-1

[Redacted] (b)(4) Trade Secret

[Redacted]

[Redacted]

[Redacted]

[Redacted]

[Redacted]

[Redacted]

[Redacted]

[Redacted]

[Redacted]

[Redacted]

# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

---

## Table of Tables

(b)(4) Trade Secret Process

[Redacted content]

# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

---

## Table of Figures

[REDACTED]	(
[REDACTED]	b
[REDACTED]	)
[REDACTED]	4
[REDACTED]	)
[REDACTED]	T
[REDACTED]	r
[REDACTED]	a
[REDACTED]	d
[REDACTED]	e



# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

---

## 1 Medical Device User Fee Cover Sheet (Form FDA 3601)

Payment Identification Number: (b)(4) Trade Secret

The Medical Device User Fee Cover Sheet (Form FDA 3601) begins on the next page.

Form Approved OMB No. 0910-0511 Expiration Date April 30, 2016. See Instructions for OMB Statement.

DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION <b>MEDICAL DEVICE USER FEE COVER SHEET</b>		PAYMENT IDENTIFICATION NUMBER: (b)(4) Trade Write the Payment Identification number on your check.	
A completed cover sheet must accompany each original application or supplement subject to fees. If payment is sent by U.S. mail or courier, please include a copy of this completed form with payment. Payment and mailing instructions can be found at: <a href="http://www.fda.gov/oc/mdufma/cover sheet.html">http://www.fda.gov/oc/mdufma/cover sheet.html</a>			
1. COMPANY NAME AND ADDRESS (include name, street address, city state, country, and post office code)  STRAUMANN USA 60 MINUTEMAN ROAD ANDOVER MA 01810 US  1.1 EMPLOYER IDENTIFICATION NUMBER (EIN) <b>(b)(4)</b>		2. CONTACT NAME Jennifer Jackson  2.1 E-MAIL ADDRESS jennifer.jackson@straumann.com  2.2 TELEPHONE NUMBER (include Area code) 978-747-2509  2.3 FACSIMILE (FAX) NUMBER (Include Area code) 978-747-0023	
3. TYPE OF PREMARKET APPLICATION (Select one of the following in each column; if you are unsure, please refer to the application descriptions at the following web site: <a href="http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm345263.htm">http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm345263.htm</a> <u>Select an application type:</u>			
<input checked="" type="checkbox"/> Premarket notification(510(k)); except for third party <input type="checkbox"/> 513(g) Request for Information <input type="checkbox"/> Biologics License Application (BLA) <input type="checkbox"/> Premarket Approval Application (PMA) <input type="checkbox"/> Modular PMA <input type="checkbox"/> Product Development Protocol (PDP) <input type="checkbox"/> Premarket Report (PMR) <input type="checkbox"/> 30-Day Notice		3.1 Select a center <input checked="" type="checkbox"/> CDRH <input type="checkbox"/> CBER  <u>3.2 Select one of the types below</u> <input checked="" type="checkbox"/> Original Application <u>Supplement Types:</u> <input type="checkbox"/> Efficacy (BLA) <input type="checkbox"/> Panel Track (PMA, PMR, PDP) <input type="checkbox"/> Real-Time (PMA, PMR, PDP) <input type="checkbox"/> 180-day (PMA, PMR, PDP)	
4. ARE YOU A SMALL BUSINESS? (See the instructions for more information on determining this status) <input type="checkbox"/> YES, I meet the small business criteria and have submitted the required qualifying documents to FDA <input checked="" type="checkbox"/> NO, I am not a small business 4.1 If Yes, please enter your Small Business Decision Number:			
5. FDA WILL NOT ACCEPT YOUR SUBMISSION IF YOUR COMPANY HAS NOT PAID AN ESTABLISHMENT REGISTRATION FEE THAT IS DUE TO FDA. HAS YOUR COMPANY PAID ALL ESTABLISHMENT REGISTRATION FEES THAT ARE DUE TO FDA? <input checked="" type="checkbox"/> YES (All of our establishments have registered and paid the fee, or this is our first device, and we will register and pay the fee within 30 days of FDA's approval/clearance of this device.) <input type="checkbox"/> NO (If "NO," FDA will not accept your submission until you have paid all fees due to FDA. This submission will not be processed; see <a href="http://www.fda.gov/cdrh/mdufma">http://www.fda.gov/cdrh/mdufma</a> for additional information)			
6. IS THIS PREMARKET APPLICATION COVERED BY ANY OF THE FOLLOWING USER FEE EXCEPTIONS? IF SO, CHECK THE APPLICABLE EXCEPTION. <input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates <input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population <input type="checkbox"/> This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only <input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially			
7. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION FOR WHICH FEES WERE WAIVED DUE TO SOLE USE IN A PEDIATRIC POPULATION THAT NOW PROPOSES CONDITION OF USE FOR ANY ADULT POPULATION? (If so, the application is subject to the fee that applies for an original premarket approval application (PMA). <input type="checkbox"/> YES <input checked="" type="checkbox"/> NO			
PAPERWORK REDUCTION ACT STATEMENT Public reporting burden for this collection of information is estimated to average 18 minutes per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to the address below.  Department of Health and Human Services, Food and Drug Administration, Office of Chief Information Officer, 1350 Piccard Drive, 4th Floor Rockville, MD 20850 [Please do NOT return this form to the above address, except as it pertains to comments on the burden estimate.]			
8. USER FEE PAYMENT AMOUNT SUBMITTED FOR THIS PREMARKET APPLICATION <b>(b)(4)</b>		03-Jul-2013	

Form FDA 3601 (01/2007)

["Close Window"](#) [Print Cover sheet](#)

# **Traditional 510(k) Submission**

Straumann® Variobase™ Abutments

---

## **2 CDRH Premarket Review Submission Cover Sheet**

The CDRH Premarket Review Submission Cover Sheet begins on the next page.

**CDRH PREMARKET REVIEW SUBMISSION COVER SHEET**

Date of Submission 07-16-2013	User Fee Payment ID Number <b>(b)(4) Trade</b>	FDA Submission Document Number (if known)
----------------------------------	---	---

**SECTION A TYPE OF SUBMISSION**

<b>PMA</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Premarket Report <input type="checkbox"/> Modular Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment <input type="checkbox"/> Licensing Agreement	<b>PMA &amp; HDE Supplement</b> <input type="checkbox"/> Regular (180 day) <input type="checkbox"/> Special <input type="checkbox"/> Panel Track (PMA Only) <input type="checkbox"/> 30-day Supplement <input type="checkbox"/> 30-day Notice <input type="checkbox"/> 135-day Supplement <input type="checkbox"/> Real-time Review <input type="checkbox"/> Amendment to PMA & HDE Supplement <input type="checkbox"/> Other	<b>PDP</b> <input type="checkbox"/> Original PDP <input type="checkbox"/> Notice of Completion <input type="checkbox"/> Amendment to PDP	<b>510(k)</b> <input checked="" type="checkbox"/> Original Submission: <input checked="" type="checkbox"/> Traditional <input type="checkbox"/> Special <input type="checkbox"/> Abbreviated (Complete section I, Page 5) <input type="checkbox"/> Additional Information <input type="checkbox"/> Third Party	<b>Request for Feedback</b> <input type="checkbox"/> Pre-Submission <input type="checkbox"/> Informational Meeting <input type="checkbox"/> Submission Issue Meeting <input type="checkbox"/> Day 100 Meeting <input type="checkbox"/> Agreement Meeting <input type="checkbox"/> Determination Meeting <input type="checkbox"/> Study Risk Determination <input type="checkbox"/> Other (specify):
<b>IDE</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement	<b>Humanitarian Device Exemption (HDE)</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment	<b>Class II Exemption Petition</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information	<b>Evaluation of Automatic Class III Designation (De Novo)</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information	<b>Other Submission</b> <input type="checkbox"/> 513(g) <input type="checkbox"/> Other (describe submission):

Have you used or cited Standards in your submission?  Yes  No (If Yes, please complete Section I, Page 5)

**SECTION B SUBMITTER, APPLICANT OR SPONSOR**

Company / Institution Name Straumann USA, LLC		Establishment Registration Number (if known) 1222315	
Division Name (if applicable)		Phone Number (including area code) 978-747-2509	
Street Address 60 Minuteman Road		FAX Number (including area code) 978-747-0023	
City Andover	State / Province MA	ZIP/Postal Code 01810	Country USA
Contact Name Jennifer M. Jackson, MS			
Contact Title Senior Regulatory Affairs Project Manager		Contact E-mail Address jennifer.jackson@straumann.com	

**SECTION C APPLICATION CORRESPONDENT (e.g., consultant, if different from above)**

Company / Institution Name			
Division Name (if applicable)		Phone Number (including area code)	
Street Address		FAX Number (including area code)	
City	State / Province	ZIP Code	Country
Contact Name			
Contact Title		Contact E-mail Address	

**SECTION D1****REASON FOR APPLICATION - PMA, PDP, OR HDE**

<input type="checkbox"/> New Device <input type="checkbox"/> Withdrawal <input type="checkbox"/> Additional or Expanded Indications <input type="checkbox"/> Request for Extension <input type="checkbox"/> Post-approval Study Protocol <input type="checkbox"/> Request for Applicant Hold <input type="checkbox"/> Request for Removal of Applicant Hold <input type="checkbox"/> Request to Remove or Add Manufacturing Site	<input type="checkbox"/> Change in design, component, or specification: <input type="checkbox"/> Software/Hardware <input type="checkbox"/> Color Additive <input type="checkbox"/> Material <input type="checkbox"/> Specifications <input type="checkbox"/> Other ( <i>specify below</i> )	<input type="checkbox"/> Location change: <input type="checkbox"/> Manufacturer <input type="checkbox"/> Sterilizer <input type="checkbox"/> Packager
<input type="checkbox"/> Process change: <input type="checkbox"/> Manufacturing <input type="checkbox"/> Packaging <input type="checkbox"/> Sterilization <input type="checkbox"/> Other ( <i>specify below</i> )	<input type="checkbox"/> Labeling change: <input type="checkbox"/> Indications <input type="checkbox"/> Instructions <input type="checkbox"/> Performance Characteristics <input type="checkbox"/> Shelf Life <input type="checkbox"/> Trade Name <input type="checkbox"/> Other ( <i>specify below</i> )	<input type="checkbox"/> Report Submission: <input type="checkbox"/> Annual or Periodic <input type="checkbox"/> Post-approval Study <input type="checkbox"/> Adverse Reaction <input type="checkbox"/> Device Defect <input type="checkbox"/> Amendment
<input type="checkbox"/> Response to FDA correspondence:		<input type="checkbox"/> Change in Ownership <input type="checkbox"/> Change in Correspondent <input type="checkbox"/> Change of Applicant Address
<input type="checkbox"/> Other Reason ( <i>specify</i> ):		

**SECTION D2****REASON FOR APPLICATION - IDE**

<input type="checkbox"/> New Device <input type="checkbox"/> New Indication <input type="checkbox"/> Addition of Institution <input type="checkbox"/> Expansion / Extension of Study <input type="checkbox"/> IRB Certification <input type="checkbox"/> Termination of Study <input type="checkbox"/> Withdrawal of Application <input type="checkbox"/> Unanticipated Adverse Effect <input type="checkbox"/> Notification of Emergency Use <input type="checkbox"/> Compassionate Use Request <input type="checkbox"/> Treatment IDE <input type="checkbox"/> Continued Access	<input type="checkbox"/> Change in: <input type="checkbox"/> Correspondent/Applicant <input type="checkbox"/> Design/Device <input type="checkbox"/> Informed Consent <input type="checkbox"/> Manufacturer <input type="checkbox"/> Manufacturing Process <input type="checkbox"/> Protocol - Feasibility <input type="checkbox"/> Protocol - Other <input type="checkbox"/> Sponsor	<input type="checkbox"/> Response to FDA Letter Concerning: <input type="checkbox"/> Conditional Approval <input type="checkbox"/> Deemed Approved <input type="checkbox"/> Deficient Final Report <input type="checkbox"/> Deficient Progress Report <input type="checkbox"/> Deficient Investigator Report <input type="checkbox"/> Disapproval <input type="checkbox"/> Request Extension of Time to Respond to FDA  <input type="checkbox"/> Request Meeting <input type="checkbox"/> Request Hearing
<input type="checkbox"/> Report submission: <input type="checkbox"/> Current Investigator <input type="checkbox"/> Annual Progress Report <input type="checkbox"/> Site Waiver Report <input type="checkbox"/> Final		
<input type="checkbox"/> Other Reason ( <i>specify</i> ):		

**SECTION D3****REASON FOR SUBMISSION - 510(k)**

<input checked="" type="checkbox"/> New Device	<input type="checkbox"/> Additional or Expanded Indications	<input type="checkbox"/> Change in Technology
<input type="checkbox"/> Other Reason ( <i>specify</i> ):		

**SECTION E ADDITIONAL INFORMATION ON 510(K) SUBMISSIONS**

Product codes of devices to which substantial equivalence is claimed								Summary of, or statement concerning, safety and effectiveness information <input checked="" type="checkbox"/> 510 (k) summary attached <input type="checkbox"/> 510 (k) statement
1	NHA	2		3		4		
5		6		7		8		

Information on devices to which substantial equivalence is claimed (if known)

	510(k) Number		Trade or Proprietary or Model Name		Manufacturer
1	K120822	1	Straumann CARES Variobase Abutment NNC, RN, WN, NC, RC	1	Straumann
2	K111935	2	Ti-Base Abutment	2	NT-Trading GmbH & Co. KG
3		3		3	
4		4		4	
5		5		5	
6		6		6	

**SECTION F PRODUCT INFORMATION - APPLICATION TO ALL APPLICATIONS**

Common or usual name or classification name  
Endosseous dental implant abutment

	Trade or Proprietary or Model Name for This Device		Model Number
1	(b)(4) Trade Secret Process	1	(b)(4)
2	(b)(4) Trade Secret Process	2	(b)(4)
3	(b)(4) Trade Secret Process	3	(b)(4)
4	(b)(4) Trade Secret Process	4	(b)(4)
5	(b)(4) Trade Secret Process	5	(b)(4)

FDA document numbers of all prior related submissions (regardless of outcome)

1	2	3	4	5	6
7	8	9	10	11	12

Data Included in Submission  
 Laboratory Testing       Animal Trials       Human Trials

**SECTION G PRODUCT CLASSIFICATION - APPLICATION TO ALL APPLICATIONS**

Product Code NHA	C.F.R. Section (if applicable) 21 CFR 872.3630	Device Class <input type="checkbox"/> Class I <input checked="" type="checkbox"/> Class II <input type="checkbox"/> Class III <input type="checkbox"/> Unclassified
Classification Panel Dental		

Indications (from labeling)  
 The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.

**Note:** Submission of the information entered in Section H does not affect the need to submit device establishment registration.

FDA Document Number *(if known)*

**SECTION H MANUFACTURING / PACKAGING / STERILIZATION SITES RELATING TO A SUBMISSION**

<input checked="" type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete	Facility Establishment Identifier (FEI) Number	<input checked="" type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Manufacturer	<input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager / Relabeler
Company / Institution Name Institut Straumann AG		Establishment Registration Number 9613348	
Division Name <i>(if applicable)</i>		Phone Number <i>(including area code)</i> 978-747-2509	
Street Address Peter Merian-Weg 12		FAX Number <i>(including area code)</i> 978-747-0023	
City Basel		State / Province	ZIP Code CH-4052 Country Switzerland
Contact Name Jennifer M. Jackson, MS	Contact Title Senior Regulatory Affairs Project Manager	Contact E-mail Address jennifer.jackson@straumann.com	

<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete	Facility Establishment Identifier (FEI) Number	<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Manufacturer	<input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager / Relabeler
Company / Institution Name		Establishment Registration Number	
Division Name <i>(if applicable)</i>		Phone Number <i>(including area code)</i>	
Street Address		FAX Number <i>(including area code)</i>	
City		State / Province	ZIP Code Country
Contact Name	Contact Title	Contact E-mail Address	

<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete	Facility Establishment Identifier (FEI) Number	<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Manufacturer	<input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager / Relabeler
Company / Institution Name		Establishment Registration Number	
Division Name <i>(if applicable)</i>		Phone Number <i>(including area code)</i>	
Street Address		FAX Number <i>(including area code)</i>	
City		State / Province	ZIP Code Country
Contact Name	Contact Title	Contact E-mail Address	

**SECTION I**

**UTILIZATION OF STANDARDS**

**Note:** Complete this section if your application or submission cites standards or includes a "Declaration of Conformity to a Recognized Standard" statement.

	Standards No.	Standards Organization	Standards Title	Version	Date
1	(b)(4) Trade Secret Process				
2	(b)(4) Trade Secret Process				
3	(b)(4) Trade Secret Process				
4	(b)(4) Trade Secret Process				
5	(b)(4) Trade Secret Process				
6	(b)(4) Trade Secret Process				
7	(b)(4) Trade Secret Process				

**Please include any additional standards to be cited on a separate page.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 0.5 hour per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
 Food and Drug Administration  
 Office of Chief Information Officer  
 Paperwork Reduction Act (PRA) Staff  
 1350 Piccard Drive, Room 400  
 Rockville, MD 20850

*An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number.*



**SECTION I**

**UTILIZATION OF STANDARDS**

**Note:** Complete this section if your application or submission cites standards or includes a "Declaration of Conformity to a Recognized Standard" statement.

	Standards No.	Standards Organization	Standards Title	Version	Date
1	<b>(b)(4) Trade Secret Process</b>				
3					
4					
5					
6					
7					

**Please include any additional standards to be cited on a separate page.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 0.5 hour per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
 Food and Drug Administration  
 Office of Chief Information Officer  
 Paperwork Reduction Act (PRA) Staff  
 1350 Piccard Drive, Room 400  
 Rockville, MD 20850

*An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number.*

# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

---

## **3 510(k) Cover Letter**

The 510(k) Cover Letter begins on the next page.



July 16, 2013

U.S. Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center – W066-G609  
10903 New Hampshire Avenue  
Silver Spring, Maryland 20993-0002

**Subject: Traditional 510(k) Premarket Notification: Straumann® Variobase™ Abutments**

Dear Sir or Madam:

On behalf of Institut Straumann AG, Straumann USA, LLC submits this Traditional 510(k) Premarket Notification for the Straumann® Variobase™ Abutments in accordance with 21 CFR 807.81. The intended use and fundamental operating principles of the proposed devices are substantially equivalent to previously cleared devices as detailed in this premarket submission.

**Submitter:**

Straumann USA, LLC (on behalf of Institut Straumann AG)  
60 Minuteman Road  
Andover, MA 01810

**Primary Contact:**

Jennifer M. Jackson, MS  
Senior Regulatory Affairs Project Manager  
Straumann USA, LLC  
60 Minuteman Road  
Andover, MA 01810  
Telephone: 800-448-8168 x2509  
Fax: 978-747-0023

**Classification Name of Device:**

Classification Name:	Abutment, Implant, Dental, Endosseous
Device Product Code:	NHA
Product Classification:	Class II
Panel:	Dental
Regulation Number:	§872.3630



**Design and Use of the Device:**

Question	YES	NO
Is the device intended for prescription use (21 CFR 801 Subpart D)?	X	
Is the device intended for over-the counter use (21 CFR 807 Subpart C)?		X
Does the device contain components derived from a tissue or other biologic source?		X
Is the device provided sterile?		X
Is the device intended for single use?	X	
Is the device a reprocessed single use device?		X
If yes, does this device type require reprocessed validation data?		
Does the device contain a drug?		X
Does the device contain a biologic?		X
Does the device use software?		X
Does the submission include clinical information?		X
Is the device implanted?		X

The term “substantially equivalent” as used herein is intended to be a determination of substantial equivalence from an FDA-regulatory point of view under the Federal Food, Drug and Cosmetic Act and relates to the fact that the product can be marketed without premarket approval or reclassification. These products may be considered distinct from a patent point of view. The term “substantially equivalent” is not applicable to, and does not diminish, any patent claim related to this product or the technology used to manufacture the product.

Information in this premarket notification is considered proprietary or trade secret or confidential commercial information. The company requests that all such information not be disclosed pursuant to 18 U.S.C. §1905, 5 U.S.C. §552, 21 U.S.C. §331(j), and all other applicable laws and regulations.

Two copies of the Traditional 510(k) Premarket Notification are enclosed. The second copy is being provided in PDF format on CD and is an exact duplicate of the hardcopy. Further, in accordance with the Medical Device User Fee and Modernization Act of 2001 (“MDUFMA”), Straumann USA, LLC has submitted the appropriate application fees. A copy of the User Fee Cover Sheet is provided with the enclosed 510(k) Premarket Notification.

We trust that the foregoing information will be sufficient to permit FDA to make a finding of substantial equivalence for the proposed Straumann® Variobase™ Abutments to the currently marketed devices as presented in this premarket notification.

Please address any questions regarding this 510(k) Premarket Notification to the undersigned.

Sincerely,

Jennifer M. Jackson, MS  
Senior Regulatory Affairs Project Manager

Enclosures

# **Traditional 510(k) Submission**

## **Straumann® Variobase™ Abutments**

---

### **4 Indications for Use Statement**

The Indications for Use Statement associated with this 510(k) is located on the following page in the required format.

## Indications for Use

510(k) Number (if known):

Device Name: Straumann® Variobase™ Abutments

Indications for Use:

The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF  
NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

---

### 5 510(k) Summary

#### 5.1 Submitter's Contact Information

Straumann USA, LLC (on behalf of Institut Straumann AG)

60 Minuteman Road

Andover, MA 01810

Phone Number: 1-978-747-2509

Fax Number: 1-978-747-0023

Contact Person: Jennifer M. Jackson, MS

Date of Submission: 16-Jul-2013

#### 5.2 Name of the Device

Trade Name: Straumann® Variobase™ Abutments

Common Name: Dental Implant Abutment

Classification Name: Abutment, Implant, Dental, Endosseous

Regulation Number: §872.3630

#### 5.3 Predicate Device(s)

K120822 – Straumann CARES Variobase Abutment NNC, RN, WN, NC, RC

K111935 – Ti-Base Abutment (NT-Trading GmbH & Co. KG)

#### 5.4 Device Description

The Straumann® Variobase™ Abutments are pre-manufactured (stock) abutments, sometimes referred to as “bonding bases”. Straumann® Variobase™ Abutments are available to fit Straumann® dental implant platforms NNC (Narrow Neck CrossFit®), RN (Regular Neck), WN (Wide Neck), NC (Narrow CrossFit®), and RC (Regular CrossFit®).

#### 5.5 Intended Use

Straumann® Variobase™ Abutments are pre-manufactured (stock) abutments intended to be placed onto Straumann dental implants to provide support for customized prosthetic restorations such as crowns and bridges. Straumann® Variobase™ Abutments are

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

---

indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.

### 5.6 Technological Characteristics

Straumann® Variobase™ Abutments are pre-manufactured (stock) abutments made from a titanium-aluminum-niobium alloy.

### 5.7 Performance Testing

The material used in the manufacture of Straumann® Variobase™ Abutments is a titanium-aluminum-niobium alloy which meets the requirements of ISO 5832-11. Bench testing was performed to evaluate the fatigue load limits of the proposed Straumann® Variobase™ Abutments. Dynamic fatigue tests were conducted in accordance to the FDA guidance document *“Guidance for Industry and FDA Staff – Class II Special Controls Guidance Document: Root-form Endosseous Dental Implants and Endosseous Dental Abutments”*.

### 5.8 Conclusion

The documentation submitted in this premarket notification demonstrates that the Straumann® Variobase™ Abutments are substantially equivalent to the predicate devices.





# **Traditional 510(k) Submission**

## **Straumann® Variobase™ Abutments**

---

### **7 Class III Summary and Certification**

This section is not applicable as the subject devices have been determined to be Class II per 21 CFR 872.3630.

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

---

### **8 Financial Certification or Disclosure Statement**

This section is not applicable as there is no clinical data being submitted to support this premarket notification.

# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

---

## 9 Declarations of Conformity and Summary Reports

(b)(4) Trade Secret Process



# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process



# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process		
A.	Applicable recognized consensus standard:	(b)(4) Trade Secret Process
B.	Requirements met?	Yes
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change
D.	Requirements that were not applicable to the device:	N/A
E.	Deviations from each applicable standard that were met including justification:	N/A
F.	Differences exist, if any, between the tested device to be marketed:	No
G.	Test laboratory:	(b)(4) Trade Secret Process

Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional                       Special                       Abbreviated

STANDARD TITLE<sup>1</sup>

**(b)(4) Trade Secret Process**

**Please answer the following questions**

Yes      No

Is this standard recognized by FDA<sup>2</sup>? .....      

FDA Recognition number<sup>3</sup> ..... **#2-156**

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report<sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
 If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
 If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
 If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
 If yes, were deviations in accordance with the FDA supplemental information sheet (SIS)<sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
 If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
 If yes, report these exclusions in the summary report table.

Is there an FDA guidance<sup>6</sup> that is associated with this standard? .....         
 If yes, was the guidance document followed in preparation of this 510k? .....      

Title of guidance: Use of International Standard ISO-10993, 'Biological Evaluation of Medical Devices Part 1: Evaluation and Testing' (Replaces #G87-1 #8294) (blue book memo), May 1, 1995

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process

A.	Applicable recognized consensus standard:	(b)(4) Trade Secret Process
B.	Requirements met?	Yes
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change
D.	Requirements that were not applicable to the device:	N/A
E.	Deviations from each applicable standard that were met including justification:	N/A
F.	Differences exist, if any, between the tested device to be marketed:	No
G.	Test laboratory:	(b)(4) Trade Secret Process



Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION  
 Traditional       Special       Abbreviated

STANDARD TITLE<sup>1</sup>  
**(b)(4) Trade Secret Process**

**Please answer the following questions** Yes      No

Is this standard recognized by FDA<sup>2</sup>? .....      

FDA Recognition number<sup>3</sup> ..... **#2-153**

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report<sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
 If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
 If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
 If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
 If yes, were deviations in accordance with the FDA supplemental information sheet (SIS)<sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
 If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
 If yes, report these exclusions in the summary report table.

Is there an FDA guidance<sup>6</sup> that is associated with this standard? .....         
 If yes, was the guidance document followed in preparation of this 510k? .....      

Title of guidance: Use of International Standard ISO-10993, 'Biological Evaluation of Medical Devices Part 1: Evaluation and Testing' (Replaces #G87-1 #8294) (blue book memo), May 1, 1995

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process

A.	Applicable recognized consensus standard:	(b)(4) Trade Secret Process
B.	Requirements met?	Yes
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change
D.	Requirements that were not applicable to the device:	N/A
E.	Deviations from each applicable standard that were met including justification:	N/A
F.	Differences exist, if any, between the tested device to be marketed:	No
G.	Test laboratory:	(b)(4) Trade Secret Process

Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION  
 Traditional       Special       Abbreviated

STANDARD TITLE<sup>1</sup>  
(b)(4) Trade Secret Process

**Please answer the following questions** Yes      No

Is this standard recognized by FDA<sup>2</sup>? .....      

FDA Recognition number<sup>3</sup> ..... **#2-135**

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report<sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
 If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
 If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
 If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
 If yes, were deviations in accordance with the FDA supplemental information sheet (SIS)<sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
 If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
 If yes, report these exclusions in the summary report table.

Is there an FDA guidance<sup>6</sup> that is associated with this standard? .....         
 If yes, was the guidance document followed in preparation of this 510k? .....      

Title of guidance: Use of International Standard ISO-10993, 'Biological Evaluation of Medical Devices Part 1: Evaluation and Testing' (Replaces #G87-1 #8294) (blue book memo), May 1, 1995

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process

A.	Applicable recognized consensus standard:	(b)(4) Trade Secret Process
B.	Requirements met?	Yes
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change
D.	Requirements that were not applicable to the device:	N/A
E.	Deviations from each applicable standard that were met including justification:	N/A
F.	Differences exist, if any, between the tested device to be marketed:	No
G.	Test laboratory:	(b)(4) Trade Secret Process

Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION  
 Traditional       Special       Abbreviated

STANDARD TITLE<sup>1</sup>  
**(b)(4) Trade Secret Process**

**Please answer the following questions** Yes      No

Is this standard recognized by FDA<sup>2</sup>? .....      

FDA Recognition number<sup>3</sup> ..... # \_\_\_\_\_

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report<sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
 If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
 If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
 If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
 If yes, were deviations in accordance with the FDA supplemental information sheet (SIS)<sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
 If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
 If yes, report these exclusions in the summary report table.

Is there an FDA guidance<sup>6</sup> that is associated with this standard? .....         
 If yes, was the guidance document followed in preparation of this 510k? .....      

Title of guidance: Use of International Standard ISO-10993, 'Biological Evaluation of Medical Devices Part 1: Evaluation and Testing' (Replaces #G87-1 #8294) (blue book memo), May 1, 1995

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process

A.	Applicable recognized consensus standard:	(b)(4) Trade Secret Process
B.	Requirements met?	Yes
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change
D.	Requirements that were not applicable to the device:	N/A
E.	Deviations from each applicable standard that were met including justification:	N/A
F.	Differences exist, if any, between the tested device to be marketed:	No
G.	Test laboratory:	(b)(4) Trade Secret Process

Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION  
 Traditional       Special       Abbreviated

STANDARD TITLE<sup>1</sup>  
**(b)(4) Trade Secret Process**

Please answer the following questions	Yes	No
Is this standard recognized by FDA <sup>2</sup> .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
FDA Recognition number <sup>3</sup> .....	<b>#8-63</b>	
Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Is a summary report <sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
If no, complete a summary report table.		
Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Does this standard include acceptance criteria? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
If no, include the results of testing in the 510(k).		
Does this standard include more than one option or selection of tests? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, report options selected in the summary report table.		
Were there any deviations or adaptations made in the use of the standard? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) <sup>5</sup> ?	<input type="checkbox"/>	<input type="checkbox"/>
Were deviations or adaptations made beyond what is specified in the FDA SIS? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, report these deviations or adaptations in the summary report table.		
Were there any exclusions from the standard? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, report these exclusions in the summary report table.		
Is there an FDA guidance <sup>6</sup> that is associated with this standard? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, was the guidance document followed in preparation of this 510k? .....	<input type="checkbox"/>	<input type="checkbox"/>
Title of guidance: _____		

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process

		(b)(4) Trade Secret Process
A.	Applicable recognized consensus standard:	(b)(4) Trade Secret Process
B.	Requirements met?	Yes
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change
D.	Requirements that were not applicable to the device:	N/A
E.	Deviations from each applicable standard that were met including justification:	N/A
F.	Differences exist, if any, between the tested device to be marketed:	No
G.	Test laboratory:	(b)(4) Trade Secret Process



Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION  
 Traditional       Special       Abbreviated

STANDARD TITLE<sup>1</sup>  
(b)(4) Trade Secret Process

Please answer the following questions	Yes	No
Is this standard recognized by FDA <sup>2</sup> .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
FDA Recognition number <sup>3</sup> ..... #_____		
Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Is a summary report <sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
If no, complete a summary report table.		
Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Does this standard include acceptance criteria? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
If no, include the results of testing in the 510(k).		
Does this standard include more than one option or selection of tests? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, report options selected in the summary report table.		
Were there any deviations or adaptations made in the use of the standard? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) <sup>5</sup> ?	<input type="checkbox"/>	<input type="checkbox"/>
Were deviations or adaptations made beyond what is specified in the FDA SIS? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, report these deviations or adaptations in the summary report table.		
Were there any exclusions from the standard? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, report these exclusions in the summary report table.		
Is there an FDA guidance <sup>6</sup> that is associated with this standard? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, was the guidance document followed in preparation of this 510k? .....	<input type="checkbox"/>	<input type="checkbox"/>
Title of guidance: _____		

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process

A.	Applicable recognized consensus standard:	(b)(4) Trade Secret Process
B.	Requirements met?	Yes
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change
D.	Requirements that were not applicable to the device:	N/A
E.	Deviations from each applicable standard that were met including justification:	N/A
F.	Differences exist, if any, between the tested device to be marketed:	No
G.	Test laboratory:	(b)(4) Trade Secret Process

Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional                       Special                       Abbreviated

STANDARD TITLE<sup>1</sup>

(b)(4) Trade Secret Process

**Please answer the following questions**

Yes      No

Is this standard recognized by FDA<sup>2</sup>? .....      

FDA Recognition number<sup>3</sup> ..... #14-261

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report<sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
 If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
 If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
 If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
 If yes, were deviations in accordance with the FDA supplemental information sheet (SIS)<sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
 If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
 If yes, report these exclusions in the summary report table.

Is there an FDA guidance<sup>6</sup> that is associated with this standard? .....         
 If yes, was the guidance document followed in preparation of this 510k? .....      

Title of guidance: \_\_\_\_\_

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process

A.	Applicable recognized consensus standard:	(b)(4) Trade Secret Process
B.	Requirements met?	Yes
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change
D.	Requirements that were not applicable to the device:	N/A
E.	Deviations from each applicable standard that were met including justification:	N/A
F.	Differences exist, if any, between the tested device to be marketed:	No
G.	Test laboratory:	(b)(4) Trade Secret Process

(b)(4) Trade Secret Process

Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional                       Special                       Abbreviated

STANDARD TITLE<sup>1</sup>

**(b)(4) Trade Secret Process**

**Please answer the following questions**

Yes      No

Is this standard recognized by FDA<sup>2</sup>? .....      

FDA Recognition number<sup>3</sup> ..... # \_\_\_\_\_

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report<sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
 If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
 If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
 If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
 If yes, were deviations in accordance with the FDA supplemental information sheet (SIS)<sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
 If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
 If yes, report these exclusions in the summary report table.

Is there an FDA guidance<sup>6</sup> that is associated with this standard? .....         
 If yes, was the guidance document followed in preparation of this 510k? .....      

Title of guidance: \_\_\_\_\_

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

(b)(4) Trade Secret  
Process

A.	Applicable recognized consensus standard:	(b)(4) Trade Secret Process
B.	Requirements met?	Yes
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change
D.	Requirements that were not applicable to the device:	N/A
E.	Deviations from each applicable standard that were met including justification:	N/A
F.	Differences exist, if any, between the tested device to be marketed:	No
G.	Test laboratory:	(b)(4) Trade Secret Process

Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION  
 Traditional       Special       Abbreviated

STANDARD TITLE<sup>1</sup>  
**(b)(4) Trade Secret Process**

**Please answer the following questions** Yes      No

Is this standard recognized by FDA<sup>2</sup>? .....      

FDA Recognition number<sup>3</sup> ..... **#4-195**

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report<sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
 If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
 If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
 If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
 If yes, were deviations in accordance with the FDA supplemental information sheet (SIS)<sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
 If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
 If yes, report these exclusions in the summary report table.

Is there an FDA guidance<sup>6</sup> that is associated with this standard? .....         
 If yes, was the guidance document followed in preparation of this 510k? .....      

Title of guidance: Guidance for Industry and FDA Staff 'Class II Special Controls Guidance Document: Root-form Endosseous Dental Implants and Endosseous Dental Implant Abutments', May 12, 2004

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

---

### 10 Device Description

The modification proposed in this premarket notification is to the Indications for Use for the previously cleared Straumann® CARES® Variobase™ Abutments for NNC, RN, WN, NC, and RC (K120822). The modified Indications for Use would allow Straumann to market the Straumann® Variobase™ Abutment as a stand-alone component. The dental laboratory would then manufacture the respective coping and/or prosthetic restoration using a burnout coping or STL model for open CAD software.

#### 10.1 Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process





# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process



### 10.2 Basal Screws

(b)(4) Trade Secret Process



# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process



# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

---

### 10.3 Accessories

(b)(4) Trade Secret Process



### 10.4 Procedure

#### 10.4.1 Restoration Design

(b)(4) Trade Secret Process



# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

---

### 10.4.2 Dental Laboratory and Prosthetic Procedures

(b)(4) Trade Secret Process

A large black rectangular redaction box covers the entire content of section 10.4.2.

### 10.5 Conclusion

(b)(4) Trade Secret Process

A large black rectangular redaction box covers the entire content of section 10.5.

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

---

### 11 Substantial Equivalence Discussion

Within the meaning of the Medical Device Amendments Act of 1976, the proposed change to the Indications for Use for the Straumann® Variobase™ Abutments in this 510(k) premarket notification are substantially equivalent to the medical devices currently in commercial distribution listed below:

- K120822, Straumann CARES Variobase Abutment NNC, RN, WN, NC, RC
- K111935, Ti-Base Abutment (NT-Trading GmbH & Co. KG)

Straumann® Variobase™ Abutments are pre-manufactured (stock) abutments intended to be placed onto Straumann dental implants to provide support for customized prosthetic reconstructions such as crowns and bridges. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.

The NT-Trading Ti-Base Abutment is a pre-manufactured abutment supplied in two parts, the abutment and screw, for fixation onto dedicated endosseous dental implants and is intended for use as an aid in prosthetic restoration. The Ti-Base is compatible with commercially available dental CAD/CAM systems, such as 3Shape, Exocad, and Dental Wings. Such systems must be validated by the user.

The NT-Trading catalog and package insert are included in this submission in Appendices 2 and 3, respectively. The products that were cleared in premarket notification K111935 are outlined in Table 5 (the 510(k) Summary is included in Appendix 4). Specifically, the L-Serie and N-Serie abutments are compatible with implants of the Straumann Dental Implant System.

(b)(4) Trade Secret Process - Product Specs



# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process - Product Specs



## Traditional 510(k) Submission

### Straumann® Variobase™ Abutments

The table below provides a comparison matrix of the proposed and predicate devices (K120822):

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
<b>K Number</b>	<b>Straumann® Variobase™ Abutments Subject Submission</b>	<b>Straumann® CARES® Variobase™ Abutments (K120822)</b>	
<b>Indications for Use</b>	<p>The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.</p>	<p>The Straumann® CARES® Variobase™ Abutment is a two-piece dental abutment consisting of the Straumann® Variobase™ Abutment and the Straumann® CARES® Variobase™ Coping which is intended to be placed onto Straumann dental implants to provide support for prosthetic reconstruction such as crown and bridges. Straumann® CARES® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.</p> <p>The Straumann® CARES® Variobase™ Coping polycon® ae in combination with the Straumann® CARES® Variobase™ Abutment is indicated for temporary (up to 180 days) dental restoration of a Straumann dental implant.</p>	<p>The indications for use are being modified to allow Straumann to market the Straumann® Variobase™ Abutment as a stand-alone component. The dental laboratory would then manufacture the respective coping and/or prosthetic restoration using a burnout coping or STL model for open CAD software.</p>
<b>Material</b>	Titanium-Aluminum-Niobium alloy (Ti-6Al-7Nb)	Titanium-Aluminum-Niobium alloy (Ti-6Al-7Nb)	Identical
<b>Abutment Diameter</b>	(b)(4) Trade Secret P d t S	(b)(4) Trade Secret Process - Product	Identical
<b>Abutment Height</b>	(b)(4) Trade Secret P d t	(b)(4) Trade Secret Process - Product	Identical

## Traditional 510(k) Submission

### Straumann® Variobase™ Abutments

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
<b>K Number</b>	<b>Straumann® Variobase™ Abutments Subject Submission</b>	<b>Straumann® CARES® Variobase™ Abutments (K120822)</b>	
<b>Mode of Action</b>	Screw-retained or cement retained	Screw-retained or cement retained	Identical
<b>Reusable</b>	No	No	Identical

**Table 6 - Comparison Matrix: Proposed Device versus Predicate Devices (K120822)**



## Traditional 510(k) Submission

### Straumann® Variobase™ Abutments

The table below provides a comparison matrix of the proposed and predicate devices (K111935):

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
<b>K Number</b>	<b>Straumann® Variobase™ Abutments Subject Submission</b>	<b>Ti-Base Abutment (K111935)</b>	
<b>Indications for Use</b>	<p>The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.</p>	<p>Ti-Base for individual Zirconium Abutments: The devices covered by this submission are abutments which are placed into a dental implant to provide support for a prosthetic restoration.</p> <p>The Ti-Base for individual Zirconium Abutments are intended for use to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support single and multiple tooth prosthesis, in the mandible or maxilla. The prosthesis can be cement retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.</p> <p>The Ti-Base abutments are indicated for use with the following implant systems:</p> <ul style="list-style-type: none"> <li>• Nobel Biocare® Replace Select®</li> <li>• Nobel Biocare NobelActive™</li> <li>• Biomet 3i® Osseotite®</li> <li>• Biomet 3i® Osseotite® Certain®</li> <li>• Nobel Biocare Branemark®</li> <li>• Straumann® synOcta®</li> <li>• Straumann® Bone Level®</li> <li>• Zimmer® Tapered Screw-vent®</li> <li>• Astra Tech OsseoSpeed®</li> <li>• Dentsply-Friadent® Frialit®</li> </ul>	Equivalent

## Traditional 510(k) Submission

### Straumann® Variobase™ Abutments

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
K Number	Straumann® Variobase™ Abutments Subject Submission	Ti-Base Abutment (K111935)	
Material	Titanium-Aluminum-Niobium alloy (Ti-6Al-7Nb)	Titanium-Aluminum-Vanadium alloy (Ti-6Al-4V)	Equivalent
Abutment Diameter	(b)(4) Trade Secret Process - Product	(b)(4) Trade Secret Process - Product Specs	Equivalent
Abutment Height	(b)(4) Trade Secret Process - Product	(b)(4) Trade Secret Process -	Equivalent
Mode of Action	Screw-retained or cement retained	Screw-retained or cement retained	Identical
Reusable	No	No	Identical

**Table 7 - Comparison Matrix: Proposed Device versus Predicate Devices (K111935)**

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### 12 Proposed Labeling

#### 12.1 Package Label

There are no changes to the package label as a result of the proposed change in this premarket notification. To aid in the review of the submission, an example of the label is shown in Figure 1.





Figure 1 - Example label for Straumann® Variobase™ Abutment

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### 12.2 Proposed Package Insert/Instructions for Use for Straumann® Variobase™ Abutments

Instructions for use: Straumann® Variobase™ Abutments




Institut Straumann AG, Peter-Merian-Weg 12, CH-4002 Basel/Switzerland,

[www.straumann.com](http://www.straumann.com)

#### English

**CAUTION: Federal law restricts this device to sale by or on the order of a dental professional.**



#### 1. Product Description

##### Abutments

Abutments are placed into dental implants to provide support for prosthetic reconstructions such as crowns and bridges.

##### Basal Screws

Basal screws are used for the fixation of the abutment to the dental implant.

#### 2. Intended use

Prosthetic components directly connected to the endosseous dental implant are intended for use as an aid in prosthetic rehabilitations.

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

---

### 3. Indications

The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.

### 4. Contraindications

Allergies or hypersensitivity to materials used, which are indicated in the following table:

Straumann Variobase™ abutments, Screws	Titanium alloy, Ti-6Al-7Nb (titanium-aluminum-niobium or TAN).

### 5. Warnings and Precaution

Our products must be secured against aspiration when used intraorally. Failure to follow the procedures outlined in these instructions may lead to any or all of the following complications:

- Aspiration or swallowing of a component
- Breakage
- Infection

The Straumann® Variobase™ abutments are single use devices.

Place implant-borne restorations only in occlusion when the implant is completely osseointegrated.

Angled abutments should not be used in areas of high mechanical loads on small diameter implants.

Dental cement or any other material used for the attachment of prosthetic components should be processed as specified by the manufacturer.

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

Implants are only to be restored with the corresponding original abutment compatible with that specific implant.

The Straumann® Variobase™ abutments have not been evaluated for safety and compatibility in the magnetic resonance environment. The Straumann® Variobase™ abutments have not been tested for heating or migration in the magnetic resonance environment.

### 6. Compatibility information

Straumann® implants and the prosthetic components are available in a variety of configurations to meet your clinical needs. The label on each product uses abbreviations to help you identify whether a particular abutment or coping is compatible with the implant that you are restoring. The implant as well as the prosthetic component contains an identifier for the connection type, summarized in the table below.

Implant connection type	Compatible prostheses
NC (Narrow CrossFit®)	parts labeled NC
RC (Regular CrossFit®)	parts labeled RC
NNC (Narrow Neck CrossFit®)	parts labeled NNC
RN (Regular Neck)	parts labeled RN
WN (Wide Neck)	parts labeled WN

### 7. Cleaning and Disinfection

Straumann® Variobase™ abutments and components are non-sterile when delivered. Before placing the restoration in the patient's mouth, the product must be cleaned, disinfected and sterilized. Straumann recommends the following procedure for cleaning, disinfection and sterilization of abutments prior to use.

- 1) Clean rinsing under flowing water while brushing outer and inner side with adequate brushes.
- 2) The pre-treated product is to be cleaned/disinfected in an automated washer disinfector. Select the appropriate program according to the manufacturer's instructions.

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### 8. Sterilization

Material	Method	Conditions
Variobase™ Abutment, TAN Screw, TAN	Autoclave (moist heat) Displacement: gravity or fractionated vacuum	134 °C (273 °F) 5 minutes

**Please note:** User should ensure the use of the appropriate biological indicator for the sterilizer and parameters used.

**Caution:** Use devices immediately after sterilization. Do not store sterilized devices.

### 9. Procedure

#### Use and handling of the Straumann® Variobase™ abutments for the Dental Technician

##### Restoration design

Make a coping or crown following standard procedure according to the material manufacturer's instructions. When using pressing or casting techniques via wax-up, use the burn-out coping for Variobase™ abutments which supports a clean and sharp-edged finish of the screw channel and a good fit to the Straumann® Variobase™ abutment.

##### Processing

Process the coping or crown following standard procedure according to the material manufacturer's instructions. Always finalize the crown or coping prior to bonding to the Straumann® Variobase™ abutments.

##### Bonding

**Please note:** It is not necessary to sandblast the Straumann® Variobase™ abutment.

- 1) Fix the abutment to the implant analog with a screw (hand-tight).
- 2) Seal the screw channel with wax.
- 3) Apply self-adhesive dental cement on the abutment. Only suitable self-adhesive cementation systems for the material used shall be used. Follow the instructions for use of both the dental material and cement/bonding material manufacturer. (Straumann® recommends Panavia™ F2.0 resin cement by Kuraray)
- 4) Bond the coping to the abutment.

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

- 5) Immediately remove excess cement from the abutment and polish the lower margin of the coping after the cement is set.
- 6) Clean the restoration prior to sending it to the dentist.
- 7) Include this instruction for use when sending the restoration to the dentist.

### Use and handling of the Straumann® Variobase™ Abutments for the Dentist

#### Remove the restoration from the master cast or the analog.

Clean, disinfect and sterilize the device as described in sections 7 and 8 of this Instructions for Use document.

#### Placing the restoration

- a) Remove the healing cap or temporary restoration.
- b) Clean and dry the interior of the implant and the abutment thoroughly.
- c) Place the sterilized restoration into the patient's mouth.
- d) Make sure that the retentive elements of the implant abutment connection are properly aligned.
- e) Use the screw delivered with titanium base or abutment to screw the final restoration into the dental implant.

**Please note:** Always ensure that the surfaces of threads and screw heads are clean and that a new screw is used for the restoration.

- f) Straumann® abutments are fixed to the implant using the Straumann® SCS screwdriver, ratchet and torque control device. Use the respective torque according to the table below:

Device type	Tightening torque	Special considerations
Abutments (permanent)	35 Ncm	n/a
Temporary abutments	15 – 35 Ncm	Tighten only to 35 Ncm if the implant is fully osseointegrated
Components on implant analogs	Hand-tight	n/a



# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

---

### **Warning**

Torques greater than 35 Ncm may result in failure of the abutment and/or implant. Torque values less than the recommended values may result in loosening of the abutment, which may lead to abutment and/or implant failure.

### **Please note**

Once the Straumann® abutment has been secured to the implant using the indicated torque, it should not be removed.

### **10. Further Information**

For additional information about the use of Straumann® products, call Straumann's customer service department or visit [www.straumann.com](http://www.straumann.com).

For additional information, consult:

Basic information on the Straumann® Variobase™ abutment

### **11. Please note**

Practitioners must have appropriate knowledge and instruction in the handling of the Straumann product described herein ("Straumann Product") for using the Straumann Product safely and properly in accordance with these instructions for use.

The Straumann Product must be used in accordance with the instructions for use provided by the manufacturer. It is the practitioner's responsibility to use the device in accordance with these instructions for use and to determine, if the device fits to the individual patient situation.

The Straumann Product is part of an overall concept and must be used only in conjunction with the corresponding original components and instruments distributed by Institut Straumann AG, its ultimate parent company and all affiliates or subsidiaries of such parent company ("Straumann"), except if stated otherwise in these instructions for use. If use of products made by third parties is not recommended by Straumann in these instructions for use, any such use will void any warranty or other obligation, express or implied, of Straumann.

### **12. Validity**

Upon publication of these instructions for use, all previous versions are superseded.

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

---

© Institut Straumann AG, 2013. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

Panavia™ is a trademark of Kurary Co, LTD, JP.

### Availability

Some items of the Straumann® Dental Implant System are not available in all countries.



Straumann Products with the CE mark fulfill the requirements of the Medical Devices Directive 93/42 EEC /



Refer to package insert



Manufacturer



Article number



Lot Number



Do not re-use



Non-sterile

Rx only

Federal law restricts this device to sale by or on the order of a dental professional.

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

---

### 13 Sterilization and Shelf Life

#### 13.1 Sterilization

(b)(4) Trade Secret Process - Product Specs



#### 7. Cleaning and Disinfection

(b)(4) Trade Secret Process - Product Specs



#### 8. Sterilization

Material	Method	Conditions
(b)(4) Trade Secret Process - Product Specs		

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process - Product Specs



### 13.2 Shelf Life

(b)(4) Trade Secret Process - Product Specs



# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

---

## 13.3 Packaging

(b)(4) Trade Secret Process - Product Specs



# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

---

## 14 Biocompatibility

(b)(4) Trade Secret Process - Product Specs



# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process - Product Specs



# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

---

### 15 Software

Software requirements do not apply to Straumann® Variobase™ Abutments and basal screws because these are standard stock products that do not contain software.

Therefore, this section is not applicable as the proposed change does not contain and is not dependent on the use of software.



# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

---

## 16 Electromagnetic Compatibility and Electrical Safety

(b)(4) Trade Secret Process - Product Specs



# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

---

### 17 Performance Testing – Bench

(b)(4) Trade Secret Process - Testing

A large black rectangular redaction box covering the entire content of section 17.

#### 17.1 Design Control

(b)(4) Trade Secret Process - Testing

A large black rectangular redaction box covering the entire content of section 17.1.

#### 17.2 Straumann® Variobase™ Abutments Material Properties

(b)(4) Trade Secret Process - Testing

A large black rectangular redaction box covering the entire content of section 17.2.

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

---

### 17.2.1 Chemical Composition TAN

(b)(4) Trade Secret Process - Testing

A large black rectangular redaction box covers the entire content of section 17.2.1.

### 17.2.2 Mechanical Characteristics TAN

(b)(4) Trade Secret Process - Testing

A large black rectangular redaction box covers the entire content of section 17.2.2.

### 17.3 Dynamic Fatigue Testing

(b)(4) Trade Secret Process - Testing

A large black rectangular redaction box covers the entire content of section 17.3.

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process - Testing



# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process - Testing



# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process - Testing



# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process - Testing



# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process - Testing





# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process - Testing



# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process - Testing



# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

---

(b)(4) Trade Secret Process - Testing



# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

---

### **18 Performance Testing – Animal**

This section is not applicable as no animal testing was performed in the development of the proposed devices.

# Traditional 510(k) Submission

Straumann® Variobase™ Abutments

---

## 19 Performance Testing – Clinical

This section is not applicable as clinical study results are not being submitted in this premarket notification.

# **Traditional 510(k) Submission**

Straumann® Variobase™ Abutments

---

## **Appendix 1 – Engineering Drawings**













# **Traditional 510(k) Submission**

Straumann® Variobase™ Abutments


---

## **Appendix 2 – NT-Trading Product Catalog**



# **MORE THAN ABUTMENT'S**

## CATALOGUE 2012/1



# MORE THAN ABUTMENT

## OPTIMIZED SAFETY AND QUALITY CONSTANT LOAD ACCORDING TO DIN EN ISO 14801

To analyze the behavior of the material and the design of the implants on so many of these force effects and to determine the load limits, they will be reproduced in simulations. The international standard ISO 14801 describes the test setup and the implementation of such a test of fatigue limit endosseous dental implants.

We have had constant load tests for endosseous dental implants carried out for nt-trading abutments by the National Institute of Research and Development and Measurement Technique in accordance with standard DIN EN ISO 14801. The load corresponded to 30° angle to the implant axis. In addition, we have, in conformity with DIN EN ISO 14801 had a finite element analysis test carried out.




**NATIONAL INSTITUTE OF RESEARCH  
AND DEVELOPMENT FOR MECHATRONICS AND  
MEASUREMENT TECHNIQUE**

**TESTING REPORT**  
No. 113 from 24.06.2010

**APPROVED,  
GENERAL MANAGER**  
Univ. Prof. PhD. eng. G. Ion GHEORGHE

**ISO 9001**  
**ISO 14001**  
**ISO 45001**

**UNI-NOR**  
F.POS - 8.2.A - 2



# CONTENTS

## IMPLANT SYSTEMS

<b>E-Serie</b> is compatible to Nobel Biocare Replace Select®	<b>4</b>
<b>F-Serie</b> is compatible to Nobel Biocare™ Nobel Active™	<b>9</b>
<b>H-Serie</b> is compatible to Biomet 3i Certain®	<b>11</b>
<b>I-Serie</b> is compatible to Biomet 3i Osseotite®	<b>13</b>
<b>K-Serie</b> is compatible to Nobel Biocare Brånemark®	<b>14</b>
<b>L-Serie</b> is compatible to Straumann Bone Level®	<b>15</b>
<b>N-Serie</b> is compatible to Straumann SynOcta®	<b>19</b>
<b>R-Serie</b> is compatible to Zimmer Tapered Screw-Vent®	<b>24</b>
<b>S-Serie</b> is compatible to Astra Tech OsseoSpeed®	<b>28</b>
<b>T-Serie</b> is compatible to Dentsply-Friadent Frialit/ Xive®	<b>31</b>
<b>Prosthetic tools</b>	<b>35</b>
<b>Dentokeep</b>	<b>36</b>
<b>Scan equipment</b>	<b>37</b>
<b>nt-milling tools and grinding tools</b>	<b>38</b>
<b>Sinter components</b>	<b>39</b>
<b>Tightening torques, Conditions of warranty</b>	<b>40</b>
<b>Notes</b>	<b>42</b>
<b>NT-Preform</b>	<b>43</b>

nt-trading GmbH & Co. KG meets the requirements of  
ISO 13485:2007 and the EU Directive 93/42/EEC  
Health Canada Recognized Registrar CMDCAS.



# E-SERIE

COMPATIBLE TO NOBEL BIOCARE REPLACE SELECT®



Platform  
3,5 mm NP

Platform  
4,3 mm RP

Platform  
5,0 mm WP

Platform  
6,0 mm

**Titanium Base**  
for individual milled  
Zirconium Abutment  
incl. screw  
Titan Grade 5



E 800

E 810

E 820

E 830

**Scan Body**  
for Titanium Base  
PEEK



E 00 W

E 10 W

E 20 W

E 30 W

**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK



E 9.3D3.500

E 9.3D4.300

E 9.3D5.000

E 9.3D6.000

**Lab Analog**  
Stainless Steel



E 50

E 51

E 52

E 53

**Abutment Screw**  
Titan Grade 5  
Recommended tightening  
torques 35 Ncm



E 60

E 61

E 61

E 61

# E-SERIE

COMPATIBLE TO NOBEL BIOCARE REPLACE SELECT®



Platform  
4,3 mm RP

Platform  
5,0 mm WP

**2-CONnect-Base Set**  
incl. 2-CONnect Cap and  
Cap-Screw  
Titan Grade 5



E 810 S

E 820 S

**2-CONnect-Base**  
Titan Grade 5  
Recommended tightening  
torques 35 Ncm



E 810 M

E 820 M

**2-CONnect Cap**  
Titan Grade 5



E 810 F

E 820 F

**Cap-Screw**  
Titan Grade 5  
Recommended tightening  
torques 30 Ncm



N 60

N 60

**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK



E 9.3D4.300

E 9.3D5.000

**Lab Analog**  
Stainless Steel



E 51

E 52



# E-SERIE

COMPATIBLE TO NOBEL BIOCARE REPLACE SELECT®



Platform  
3,5 mm NP

Platform  
4,3 mm RP

Platform  
5,0 mm WP



## Straight Abutment

incl. Screw  
GH 1,0 mm

E 100-1

E 110-1

E 120-1

## Straight Abutment

incl. Screw  
GH 2,5 mm

E 100

E 110

E 120

# E-SERIE

COMPATIBLE TO NOBEL BIOCARE REPLACE SELECT®



Platform  
3,5 mm NP

Platform  
4,3 mm RP

Platform  
5,0 mm WP



## Angled Abutment 16°

angled over surface  
GH 1,0 mm  
incl. Screw

E 200-1-1

E 210-1-1

E 220-1-1



## Angled Abutment 16°

angled over edge  
GH 1,0 mm  
incl. Screw

E 200-2-1

E 210-2-1

E 220-2-1



## Angled Abutment 16°

angled over surface  
GH 2,5 mm  
incl. Screw

E 200-1

E 210-1

E 220-1



## Angled Abutment 16°

angled over edge  
GH 2,5 mm  
incl. Screw

E 200-2

E 210-2

E 220-2



# E-SERIE

COMPATIBLE TO NOBEL BIOCARE REPLACE SELECT®



Platform  
3,5 mm NP

Platform  
4,3 mm RP

Platform  
5,0 mm WP



**HSL Abutment**  
rotation indexed  
incl. Screw

E 11.CA3.500

E 11.CA4.300

E 11.CA5.000

# F-SERIE

COMPATIBLE TO NOBEL BIOCARE™ NOBEL ACTIVE™



Platform  
3,5 mm NP

Platform  
4,3 mm / 5,0 mm RP



**Titanium Base**  
for individual milled  
Zirconium Abutment  
incl. screw  
Titan Grade 5

F 800

F 810



**Scan Body**  
for Titanium Base  
PEEK

F 00 W


F 10 W



**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK

F 9.3D3.500

F 9.3D4.350



**Lab Analog**  
Stainless Steel

F 50

F 51



**Abutment Screw**  
Titan Grade 5  
Recommended tightening  
torques 35 Ncm

F 60

F 61

**Implant Pic-Up**  
incl. screw

E TR-NP013.5

E TR-RP024.3

E TR-WP035.0

# F-SERIE

KOMPATIBEL ZU NOBEL BIOCARE™ NOBEL ACTIVE™



Platform  
3,5 mm NP



Platform  
4,3 mm / 5,0 mm RP



**2-COONnect-Base Set**  
incl. 2-COONnect Cap and  
Cap-Screw  
Titan Grade 5

F 800 S


F 810 S



**2-COONnect-Base**  
Titan Grade 5  
Recommended tightening  
torques 35 Ncm

F 800 M

F 810 M



**2-COONnect Cap**  
Titan Grade 5

F 800 F

F 810 F

**Cap-Screw**  
Titan Grade 5  
Recommended tightening  
torques 30 Ncm

N 60

N 60

**Scan Body 3D Guide**  
for Titanium Base +  
2-COONnect  
PEEK

F 9.3D3.500

F 9.3D4.350



**Lab Analog**  
Stainless Steel

F 50

F 51

# H-SERIE

KOMPATIBEL ZU BIOMET 3I CERTAIN®



Platform  
3,4 mm



Platform  
4,1 mm



Platform  
5,0 mm



**Titanium Base**  
for individual milled  
Zirconium Abutment  
incl. screw  
Titan Grade 5

H 800

H 810

H 820



**Scan Body**  
for Titanium Base  
PEEK

H 00 W

H 10 W

H 10 W



**Scan Body 3D Guide**  
for Titanium Base +  
2-COONnect  
PEEK

H 9.3D3.400

H 9.3D4.150

H 9.3D4.150



**Lab Analog**  
Stainless Steel

H 50

H 51

H 52



**Abutment Screw**  
Titan Grade 5  
Recommended tightening  
torques 20 Ncm

H 60

H 60

H 60



# H-SERIE

COMPATIBLE TO BIOMET 3I CERTAIN®



Platform  
3,4 mm



Platform  
4,1 mm



**2-COConnect-Base Set**  
incl. 2-COConnect Cap and  
Cap-Screw  
Titan Grade 5

H 800 S

H 810 S



**2-COConnect-Base**  
Titan Grade 5  
Recommended tightening  
torques 20 Ncm

H 800 M

H 810 M



**2-COConnect Cap**  
Titan Grade 5

H 800 F

H 810 F



**Cap-Screw**  
Titan Grade 5  
Recommended tightening  
torques 15 Ncm

KS 60



**Scan Body 3D Guide**  
for Titanium Base +  
2-COConnect  
PEEK

H 9.3D3.400

H 9.3D4.150



**Lab Analog**  
Stainless Steel

H 50

H 51

# I-SERIE

COMPATIBLE TO BIOMET 3I OSSEOTITE®



Platform  
3,4 mm



Platform  
4,1 mm



Platform  
5,0 mm



**Titanium Base**  
for individual milled  
Zirconium Abutment  
incl. screw  
Titan Grade 5

I 800

I 810

I 820



**Scan Body**  
for Titanium Base  
PEEK

I 100 W

I 10 W

I 10 W



**Scan Body 3D Guide**  
for Titanium Base +  
2-COConnect  
PEEK

I 9.3D3.400

I 9.3D4.150

I 9.3D4.150



**Lab Analog**  
Stainless Steel

I 50

I 51

I 52



**Straight Abutment**  
incl. Screw  
GH 2,5 mm

I 110



**Abutment Screw**  
Titan Grade 5  
Recommended tightening  
torques 35 Ncm

I 61

I 61

I 61

# K-SERIE

COMPATIBLE TO NOBEL BIOACARE BRÄNEMARK®



Platform  
3,5 mm

Platform  
4,1 mm

Platform  
5,1 mm

**Titanium Base**  
for individual milled  
Zirconium Abutment  
incl. screw  
Titan Grade 5



K 800

K 810

K 820

**Scan Body**  
for Titanium Base  
PEEK



K 00W

K 10W

K 20W

**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK



K 9.3D3.500

K 9.3D4.100

K 9.3D5.100

**Lab Analog**  
Stainless Steel



K 50

K 51

K 52

**Straight Abutment**  
incl. Screw  
GH 2,5 mm



K 110

**Abutment Screw**  
Titan Grade 5  
Recommended tightening  
torques 35 Ncm



K 60

K 61

K 62

# L-SERIE


COMPATIBLE TO STRAUMANN BONE LEVEL®



Platform  
3,3 mm NC

Platform  
4,1 mm / 4,8 mm RC


**Titanium Base**  
for individual milled  
Zirconium Abutment  
incl. screw  
Titan Grade 5



L 800

L 810

**Scan Body**  
for Titanium Base  
PEEK



L 00W

L 10W

**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK



L 9.3D3.300

L 9.3D4.148

**Lab Analog**  
Stainless Steel



L 50

L 51

**Abutment Screw**  
Titan Grade 5  
Recommended tightening  
torques 35 Ncm



L 60

L 60



# L-SERIE

COMPATIBLE TO STRAJMANN BONE LEVEL®



Platform  
3,3 mm NC



Platform  
4,1 mm / 4,8 mm RC



**2-CONNECT-Base Set**  
incl. 2-CONNECT Cap and  
Cap-Screw  
Titan Grade 5

L 800 S


L 810 S



**2-CONNECT-Base**  
Titan Grade 5  
Recommended tightening  
torques 35 Ncm

L 800 M

L 810 M



**2-CONNECT Cap**  
Titan Grade 5

L 800 F

L 810 F



**Cap-Screw**  
Titan Grade 5  
Recommended tightening  
torques 30 Ncm

N 60

N 60



**Scan Body 3D Guide**  
for Titanium Base +  
2-CONNECT  
PEEK

L 9.3D3.300

L 9.3D4.148



**Lab Analog**  
Stainless Steel

L 50

L 51

# L-SERIE

COMPATIBLE TO STRAJMANN BONE LEVEL®



Platform  
3,3 mm NC



Platform  
4,1 mm / 4,8 mm RC



**Straight Abutment**  
incl. Screw  
GH 2,5 mm

L 100

L 110



**Straight Abutment**  
incl. Screw  
GH 3,0 mm

L 100-3

L 110-3



**Angled Abutment 18°**  
angled over surface  
GH 1,5 mm  
incl. Screw

L 200-1

L 210-1



**Angled Abutment 18°**  
angled over edge  
GH 1,5 mm  
incl. Screw

L 200-2

L 210-2

# L-SERIE

COMPATIBLE TO STRAUMANN BONE LEVEL®



Platform  
3,3 mm NC



Platform  
4,1 mm / 4,8 mm RC

**Angled Abutment 18°**  
angled over surface  
GH 3,0 mm  
incl. Screw



L 200-1-3



**Angled Abutment 18°**  
angled over edge  
GH 3,0 mm  
incl. Screw




L 200-2-3



**HSL Abutment**  
rotation indexed  
incl. Screw



L 11.CA3.300



L 11.CA4.148

**Implant Pic-Up**  
incl. screw



L TR-NC013.3



L TR-RC024.1

# N-SERIE

COMPATIBLE TO STRAUMANN SYNOCTA®



Platform  
3,5 mm NN



Platform  
4,8 mm RN



Platform  
6,5 mm WN


**Titanium Base**  
for individual milled  
Zirconium Abutment  
incl. screw  
Titan Grade 5



N 800



N 810



N 820

**Scan Body**  
for Titanium Base  
PEEK



N 00W



N 10W



N 20W

**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK



N 9.3D3.500



N 9.3D4.800



N 9.3D6.500

**Lab Analog**  
Stainless Steel



N 50



N 51



N 52

**Abutment Screw**  
Titan Grade 5  
Recommended tightening  
torques 35 Ncm



N 60



N 62



N 62



# N-SERIE

COMPATIBLE TO STRAUMANN SYNOCTA®




Platform  
4,8 mm RN



Platform  
6,5 mm WN

## 2-CONnect-Base Set

incl. 2-CONnect Cap and  
Cap-Screw  
Titan Grade 5



N 820 S

## 2-CONnect-Base

Titan Grade 5  
Recommended tightening  
torques 35 Ncm



N 810 M


N 820 M

## 2-CONnect Cap

Titan Grade 5

## Cap-Screw

Titan Grade 5  
Recommended tightening  
torques 30 Ncm



N 810

N 820 F

N 60

N 60

## Scan Body 3D Guide

for Titanium Base +  
2-CONnect  
PEEK



N 9.3D4.800

N 9.3D6.500

## Lab Analog

Stainless Steel



N 51

N 52

# N-SERIE

COMPATIBLE TO STRAUMANN SYNOCTA®



Platform  
4,8 mm RN



Platform  
4,8 mm RN



Platform  
6,5 mm WN

## Straight Abutment

incl. Screw




N 110

N 110 L

N 120

## Angled Abutment 16°

angled over surface  
incl. Screw




N 210-1

N 210-1

N 220-1

## Angled Abutment 16°

angled over edge  
incl. Screw




N 210-2

N 210-2

N 220-2

## Abutment Screw

Titan Grade 5  
Recom. tight. torques 35 Ncm



N 61

N 61

N 61



Platform 4,8 mm RN

Total height  
4,0 mm

Total height  
5,5 mm

Total height  
7,0 mm

## RN-Massiv Abutment

only for Dentist



N 110-40

N 110-55

N 110-70



Platform 6,5 mm WN

Total height  
5,5 mm

Total height  
5,5 mm

## WN-Massiv Abutment

only for Dentist



N 120-55

# N-SERIE

COMPATIBLE TO STRAUMANN SYNOCTA®



Platform  
3,5 mm NN

Platform  
4,8 mm RN

Platform  
6,5 mm WN

# N-SERIE

COMPATIBLE TO STRAUMANN SYNOCTA®



Platform  
3,5 mm NN

Platform  
4,8 mm RN

Platform  
6,5 mm WN



**Straight Abutment**  
incl. Screw

N 100



**Angled Abutment 16°**  
incl. Screw

N 200



**Abutment Screw**

N 60



**HSL Abutment**  
rotation indexed  
incl. Screw

N 300



**HSL Abutment**  
rotating  
incl. Screw

N 300 R



**Abutment Screw**  
Titan Grade 5  
Recommended tightening  
torques 35 Ncm

N 60



**Angled Abutment 21°**  
angled over surface  
incl. Screw

N 210-1-21



**Angled Abutment 21°**  
angled over edge  
incl. Screw

N 210-2-21



**WN-Massiv Abutment**  
Height 4,0 mm

N 120-40

N 120-40



**Implant Pic-Up**  
incl. screw

N TR-NN013.5

N TR-FND24.8

N TRAWN036.5



# R-SERIE

COMPATIBLE TO ZIMMER TAPERED SCREW-VENT®



Platform  
3,5 mm

Platform  
4,5 mm

Platform  
5,7 mm

**Titanium Base**  
for individual milled  
Zirconium Abutment  
incl. screw  
Titan Grade 5



R 800

R 810

R 820

**Scan Body**  
for Titanium Base  
PEEK



R 00W

R 10W

R 20W

**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK



R 9.3D3.500

R 9.3D4.500

R 9.3D5.700

**Lab Analog**  
Stainless Steel



R 50

R 51

R 52

**Abutment Screw**  
Hex 0,50" (1,26 mm)  
Recommended tightening  
torques 30 Ncm



R 60

R 60

R 60

# R-SERIE


COMPATIBLE TO ZIMMER TAPERED SCREW-VENT®



Platform  
3,5 mm

Platform  
4,5 mm

**2-CONnect-Base Set**  
incl. 2-CONnect Cap and  
Cap-Screw  
Titan Grade 5



R 800 S

R 810 S

**2-CONnect-Base**  
Titan Grade 5  
Recommended tightening  
torques 30 Ncm



R 800 M

R 810 M

**2-CONnect Cap**  
Titan Grade 5



R 800 F

R 810 F

**Cap-Screw**  
Titan Grade 5  
Recommended tightening  
torques 25 Ncm



KS 60

KS 60

**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK



R 9.3D3.500

R 9.3D4.500

**Lab Analog**  
Stainless Steel



R 50

R 51

# R-SERIE

COMPATIBLE TO ZIMMER TAPERED SCREW-VENT®



Platform  
3,5 mm

Platform  
4,5 mm

Platform  
5,7 mm

## Straight Abutment

incl. Screw  
GH 2,5 mm



## Angled Abutment 16°

angled over surface  
incl. Screw  
GH 2,5 mm



## Angled Abutment 16°

angled over surface  
incl. Screw  
GH 2,5 mm



## Massiv Abutment

Titan Grade 5



## HSL Abutment

rotation indexed  
incl. Screw



## Abutment Screw

Hex 0,50" (1,26 mm)  
Recommended tightening  
torques 30 Ncm



# R-SERIE

COMPATIBLE TO ZIMMER TAPERED SCREW-VENT®



Platform  
3,5 mm

Platform  
4,5 mm

## Implant Pic-Up

incl. screw





# S-SERIE

COMPATIBLE TO ASTRA TECH OSSEOSPEED®



Platform 3,5 mm / 4,0 mm 4,5 mm / 5,0 mm



**Titanium Base**  
for individual milled  
Zirconium Abutment  
incl. screw  
Titan Grade 5

S 800

S 820



**Scan Body**  
for Titanium Base  
PEEK

S 00 W

S 20 W



**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK

S 9.3D3.540

S 9.3D4.550



**Lab Analog**  
Stainless Steel

S 50

S 52



**Straight Abutment**  
incl. Screw  
GH 1,5 mm

S 100



**Abutment Screw**  
Titan Grade 5  
Recommended tightening  
torques 25 Ncm

S 60

S 61

# S-SERIE

COMPATIBLE TO ASTRA TECH OSSEOSPEED®



Platform 3,5 mm / 4,0 mm 4,5 mm / 5,0 mm



**2-CONnect-Base Set**  
incl. 2-CONnect Cap and  
Cap-Screw  
Titan Grade 5

S 800 S

S 820 S



**2-CONnect-Base**  
Titan Grade 5  
Recommended tightening  
torques 25 Ncm

S 800 M

S 820 M



**2-CONnect Cap**  
Titan Grade 5

S 800 F

S 820 F



**Cap-Screw**  
Titan Grade 5  
Recommended tightening  
torques 20 Ncm

N 60

N 60



**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK

S 9.3D3.540

S 9.3D4.550



**Lab Analog**  
Stainless Steel

S 50

S 52

## S-SERIE

COMPATIBLE TO ASTRA TECH OSSEOSPEED®



Platform 3,5 mm / 4,0 mm  
Platform 3,5 mm / 4,0 mm  
Platform 4,5 mm / 5,0 mm

<b>Straight Abutment</b> incl. Screw GH 1,5 mm			
	S 110	S 110	S 120

<b>Angled Abutment 16°</b> incl. Screw GH 1,5 mm			
	S 200	S 210	S 220

<b>HSL Abutment</b> rotation indexed incl. Screw		
	S 11,CA3,540	S 11,CA4,550

<b>Implant Pic-Up</b> incl. screw		
	STR-00013.5	STR-00024.5

## T-SERIE

COMPATIBLE TO DENTSPLY-FRIADENT FRIALIT/XIVE®



Platform 3,4 mm  
Platform 3,8 mm  
Platform 4,5 mm  
Platform 5,5 mm

<b>Titanium Base</b> for individual milled Zirconium Abutment incl. screw Titan Grade 5				
	T 800	T 805	T 810	T 820

<b>Scan Body</b> for Titanium Base PEEK				
	T 00W	T 05W	T 10W	T 10W

<b>Scan Body 3D Guide</b> for Titanium Base + 2-CONnect PEEK				
	T 9.3D3,400	T 9.3D3,800	T 9.3D4,555	T 9.3D4,555


<b>Lab Analog</b> Stainless Steel				
	T 50	T 55	T 51	T 52

<b>Abutment Screw</b> Titan Grade 5 Recommended tightening torques 25 Ncm				
	T 60	T 60	T 60	T 60



# T-SERIE

COMPATIBLE TO DENTSPLY-FRIADENT FRIALIT/XIVE®



Platform  
3,4 mm



Platform  
3,8 mm



Platform  
4,5 mm

## 2-COINnect-Base Set

incl. 2-COINnect Cap and  
Cap-Screw  
Titan Grade 5



## 2-COINnect-Base

Titan Grade 5  
Recommended tightening  
torques 25 Ncm



## 2-COINnect Cap

Titan Grade 5



## Cap-Screw

Titan Grade 5  
Recommended tightening  
torques 15 Ncm



## Scan Body 3D Guide

for Titanium Base +  
2-COINnect  
PEEK



## Lab Analog

Stainless Steel



# T-SERIE

COMPATIBLE TO DENTSPLY-FRIADENT FRIALIT/XIVE®



Platform  
4,5 mm



Platform  
3,8 mm

## Angled Abutment 16°

angled over edge  
GH 1,0 mm



## Angled Abutment 16°

angled over surface  
GH 1,0 mm



## Angled Abutment 16°

angled over surface  
GH 2,5 mm



## Straight Abutment

incl. Screw  
GH 1,0 mm



## Straight Abutment

incl. Screw  
GH 2,5 mm



## Abutment Screw

Hex 0,50" (1,26 mm)  
Recommended tightening  
torques 25 Ncm



# T-SERIE

COMPATIBLE TO DENTSPLY-FRIADENT FRIALIT/XIVE®



Platform  
3,4 mm

Platform  
3,8 mm

Platform  
4,5 mm



**Implant Pic-Up**  
incl. screw

T TR-00013.4

T TR-00023.8

T TR-00034.5



**HSL Abutment**  
notation indexed  
incl. Screw

T 11.CA3.400

T 11.CA3.800

T 11.CA4.500



# PROSTHETIC TOOLS



**Torque Ratchet**  
continuously adjustable  
to max. 40 Ncm  
ISO-Shaft Connection

W 11.000.000



**Laboratory**  
Screwdriver  
for changeable inserts  
(without ISO-Shaft)  
ISO-Shaft Connection

W 11.100.000



**Insert**  
for Laboratory Screwdriver  
and Torque Ratchet  
ISO-Shaft for S, J and R-Series  
Hex 0,50" (1,26 mm)

W 11.IFS.G10



**Insert**  
for Laboratory Screwdriver  
and Torque Ratchet  
ISO-Shaft for T- and H-Series  
Hex 1,20 mm

W 11.TH0.G20



**Insert**  
for Laboratory Screwdriver  
and Torque Ratchet  
ISO-Shaft for Ankylos®  
Hex 1,00 mm

W 11.Y00.G30



**Insert**  
for Laboratory Screwdriver  
and Torque Ratchet  
ISO-Shaft for L- and N-Series  
Torx T6

W 11.LN0.G40



**Insert**  
for Laboratory Screwdriver  
and Torque Ratchet  
ISO-Shaft for E, F- and  
K-Series UG

W 11.EFK.G50



**Insert**  
for Laboratory Screwdriver  
and Torque Ratchet ISO-Shaft  
for 2-CONNECT-Base-Primary

W 11.005.G60





# DENTOKEEP PEEK DISC

## Dentokeep

nt-trading Dentokeep Disc PEEK are blanks for technical milling manufacture of permanent dentures and prosthetic brace in the CAM process. The derived prosthetic designs are available on the remaining teeth, - can be attached to stumps and / or implant abutment and provide a basic functional and aesthetic care.

1. PEEK Disc 98.5 mm, pearl white, 14 mm and 18 mm thick
2. For CAD / CAM Technology for CAD/CAM Technology
3. Clip in prosthetics CAD/CAM
4. Class II product for permanent prosthetic restorations



## Dentokeep

PEEK Disc  
Size 14 mm

12.000.P14



## Dentokeep

PEEK Disc  
Size 18 mm

12.000.P18



# SCAN EQUIPMENT

## nt-OptiScan™ Spray

nt-OptiScan™ Spray a new product development of nt-trading GmbH & Co. KG.

## nt-OptiScan™ Spray

### Advantages:

- Optimize composition
- Perfect dosage with proven applicator
- Optimize preparation of your working model for optical scanning with modern dental scanner
- Economical application
- Water soluble
- Handy size: Volume 75 ml brutto
- Superior scanning results
- Special attractive price





# NT-MILLING AND GRINDING TOOLS

nt-trading tools are specially developed for dental CAD/CAM systems. They are up to the high mark to process dental materials and are constant liable to quality checks.

Tools coated with tax are especially qualified for process of dental high-capacity plastics, as well as for process of milling wax.

Tools coated with diamond are especially qualified for process of abrasive materials like zirconium dioxide.

<b>nt-Diaburr</b>			
Ø	0,6 mm	1,0 mm	2,5 mm
Art.-Nr:	11100	11101	11102

The milling tool coated with diamonds nt-Diaburr offers you a lifetime of at least 350 units while process of ZrO2 with the milling machine.

<b>nt-Uburr</b>			
Ø	1,0 mm	2,5 mm	
Art.-Nr:	11201	11202	

Universal milling tools are qualified for process of plastics, waxes and ZrO2 with the milling machine. Depending on choice of material the lifetime is at least 175 unites.

<b>nt-PEEKburr</b>			
Ø	1,0 mm	2,5 mm	
Art.-Nr:	11301	11302	

Special milling tools are qualified for process of PEEK materials like Dentokeep, plastics and waxes.

# SINTER COMPONENTS

## nt-pearls

Sinterization of Zr-oxide substructures.

## Sinter pearls for dental CAD/CAM technology

### Advantages:

Sinterization of Zr-oxide substructures. Full ceramics are the answer to patients demands for highly esthetic, metal-free and durable restorations Zr-oxide became a common substructure for dental restorations already quite a while ago. This material features outstanding biocompatibility and mechanical properties Zr-oxide is mechanically machined in its greenstage, taking into consideration the specific shrinkage of this material. The sinterization is optimized for the specific material and its oversized milling relative to the material properties, thus assuring the best possible fit.

Another important factor contributing to a perfect fit is the right choice of support for this to be sintered substructure. Every manufacturer has his own recommendations concerning these supports. Some suggest to sinter on special „pearls“, unfortunately particle size and quality are not always right to by example avoid getting pinched between the interproximals and thus distorting the bridge.

Other solutions favour the use of „drops“, resting on very expensive and fragile support discs. The necessary preparations for this type of support are always time consuming, sometimes even cumbersome and always generating additional costs.

A viable alternative are special hi- density and quality sinterization pearls from nt-trading, optimized for crowns and bridges and therefore avoiding the common „pinching“ and distortion problems especially known to bridges.





# INSTRUCTION FOR USE

## Indication:

For manufacturing of individual abutments on dental implants. The individual abutments can be combined with copings, crowns or suprastructures made of dental ceramics.

## Contraindication:

The Ti-Bases of each Series can only be combined with the matching implant, e.g. the E-Series shall be combined exclusively with Replace Select® Implants. They cannot be combined with implants of a different implant type or manufacturer. The diameter of the Ti-base must correspond in size to the used implant in order to prevent a peri-implant tissue irritation. The Ti-Bases are indicated for single use only. If they are used multiple times, they might damage the implants.

For fixation of the Ti-Bases on the implant, the correct torque force, recommended by the implant manufacturer, has to be considered carefully to avoid the damage of the implant-bone connection.

Ncm	Abutment
20	H-Serie
25	S-Serie T-Serie F-Serie
30	R-Serie
35	I-Serie K-Serie N-Serie E-Serie L-Serie

Mechanical treatment of the connection part of the Ti-Base will damage the correct fitting of the Ti-Base on the implant.

## Handling method:

**Ceramic abutments:**  
Milling with CAD/CAM-machines of zirconium oxide- or aluminum oxide- ceramics according to the anatomic form of a crown or coping. The ceramic copings or crowns shall be milled or polished with diamond instruments and with minimal pressure and water-cooling. The minimal thickness shall be 0.5 mm. Sharpe edges must be avoided.

## Veneering:

Copings shall be veneered with appropriate ceramics before cementing onto the Ti-Base. The instructions for use of the ceramic manufacturers have to be considered. Treatment of the Ti-Base and the ceramic abutment before cementing: Sandblasting of the contact surfaces with Al<sub>2</sub>O<sub>3</sub>, 50 µm, 2 bar and intensive cleaning of dust and grease. It is recommended to protect the connection part of the Ti-Base with an implant analog during handling.

## Cementing:

It is recommended to cement the ceramic abutment onto the Ti-Base with Panavia® F2.0 (Kuraray) with RelayXUnicem® (3M Espe) or other equivalent cements. The instructions for use of the cements shall be followed carefully. The Ti-Base shall be fixed onto an implant analog with the abutments screw. The head of the screw has to be covered with wax or resin. The mixed cement is applied onto the contact part of the Ti-Base. The abutment is pressed onto the Ti-Base. The final position is evaluated by slight rotation. The gap between abutment and the Ti-Base must be as small as possible. Remaining cement shall be removed immediately.

## Polishing:

After hardening the remaining cement shall be removed with rotating silicon instruments. The cement inside the screw channel has to be removed carefully.

## Scan Body:

### Indication:

#### Scan Body:

For the CAD/CAM scanning of the model, the Scan Body is used to indicate the position of the implant. The size of the Scan Body shall be corresponding to the original Implant system, implant diameter and Ti-Base Series. The diameter of the Scan Body prevents the rotation of the ceramic abutment. The Scan Body is fixed on the implant analogue with the abutment screw. After correct positioning, there is no gap visible between implant and Scan Body. Rotation of the Scan Body is impossible

## Tightening torques

Ncm	2-CONnect-M-Abutment	
20	H-Serie	
25	S-Serie	T-Serie
30	R-Serie	
35	E-Serie	F-Serie L-Serie N-Serie

Ncm	2-CONnect Cap-Screw N 60	
20	S-Serie	T 805, T 810
30	E-Serie	F-Serie L-Serie N-Serie

Ncm	2-CONnect Cap-Screw KS 60	
15	H 800, H 810	T 800 R 800 R 810

## Conditions of warranty

Within our general terms of sale and warranty we ensure the perfect quality of our products. Due to our high production standards can offer you a '10-years' warranty on our prosthetic components.

We offer you a warranty on our mtrading components according to the conditions of warranty.

The warranty includes all material and manufacturing defects which may occur within the '10-years' time of warranty. We only give warranty for our contracting/purchasing partners (dentists, dental hospitals, laboratory). Any other persons besides those mentioned cannot lay claims to the warranty. It is not possible to assign the warranty claims.





**The new address  
from February 2013:**

nt-trading GmbH & Co. KG  
Nördliche Uferstraße 8  
76189 Karlsruhe  
Germany

# **Traditional 510(k) Submission**

Straumann® Variobase™ Abutments

---

## **Appendix 3 – NT-Trading Ti-Base Package Insert**

## Instruction for use

### Ti- Base-, for Individual Abutments, Scan Base and Scan Body E-, I-, K-, N-, R-, S-, T-, H-, L-, and F- Series

CE 0123

#### INDICATION:

The Ti-Base, is intended for use to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support single or multiple tooth prosthesis, in the mandible or maxilla. The prosthesis can be cement retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.

The Ti-Base are for manufacturing of individual abutments on dental implants. The individual abutments can be combined with restorations such as, crowns or superstructures made of dental ceramics.

The Ti-Bases of the E-Series are indicated for Replace select<sup>®</sup> implants, manufactured by Nobel Biocare<sup>®</sup>.

The Ti-Base of the I-Series are indicated for Osseotite<sup>®</sup> implants manufactured by Biomet 3i<sup>®</sup>.

The Ti-Bases of the K-Series are indicated for Brånemark<sup>®</sup> implants, manufactured by Nobel Biocare<sup>®</sup>.

The Ti-Bases of the N-Series are indicated for synOcta<sup>®</sup> implants, manufactured by Straumann<sup>®</sup>.

The Ti-Bases of the R-Series are indicated for Tapered Screw-vent<sup>®</sup> implants, manufactured by Zimmer<sup>®</sup>.

The Ti-Bases of the S-Series are indicated for OsseoSpeed<sup>®</sup> implants, manufactured by Astra Tech<sup>®</sup>.

The Ti-Bases of the T-Series are indicated for Frialit<sup>®</sup> implants, manufactured by Dentsply-Friadent<sup>®</sup>.

The Ti-Bases of the H-Series are indicated for Biomet 3i Certain<sup>®</sup> implants, manufactured by Biomet 3i<sup>®</sup>.

The Ti-Bases of the L-Series are indicated for Bone Level implants, manufactured by Straumann<sup>®</sup>.

The Ti-Bases of the F-Series are indicated for NobelActive<sup>™</sup> implants, manufactured by Nobel Biocare<sup>®</sup>.

#### ARTICLE NO. OF TI-BASES AND ABUTMENT SCREWS

The Abutments Ti-Base E-Series is compatible with Nobel Biocare Replace Select<sup>®</sup> Implant.

Nobel Replace Select <sup>®</sup>	3.5 mm	4.3 mm	5.0 mm	6.0 mm
Ti-Base E-Series	E 800	E 810	E 820	E 830
Abutment screw	E 60	E 61	E 61	E 61

The Abutment Ti-Base I-Series is compatible with Biomet 3i Osseotite<sup>®</sup>.

Biomet 3i Osseotite <sup>®</sup>	3.4 mm	4.1 mm	5.0 mm
Ti-Base I-Series	I 800	I 810	I 820
Abutment screw	I 61	I 61	I 61

The Abutment Ti-Base K-Series is compatible with Nobel Biocare Brånemark<sup>®</sup> Implant.

Brånemark <sup>®</sup>	3.5 mm	4.1 mm	5.1 mm
Ti-Base K-Series	K 800	K 810	K 820
Abutment screw	K 60	K 61	K 62

The Abutment Ti-Base N-Series is compatible with Straumann SynOcta<sup>®</sup> Implant.

Straumann SynOcta <sup>®</sup>	3.5 mm	4.8 mm	6.5 mm
Ti-Base N-Series	N 800	N 810	N 820
Abutment screw	N 60	N 62	N 62

The Abutment Ti-Base R-Series is compatible with Zimmer Tapered Screw-vent<sup>®</sup> (Sulzer) Implant.

Zimmer (Sulzer) Tapered Screw-vent <sup>®</sup>	3.5 mm	4.5 mm	5.7 mm
Ti-Base R-Series	R 800	R 810	R 820
Abutment screw	R 60	R 60	R 60

The Abutment Ti-Base S-Series is compatible with Astra Tech OsseoSpeed<sup>®</sup> Implant.

Astra Tech Osseo-Speed <sup>®</sup>	3.5 / 4.0 mm	4.5 / 5.0 mm
Ti-Base S-Series	S 800 / 810	S 820
Abutment screw	S 60	S 61

The Abutment Ti-Base T-Series is compatible with Dentsply-Friadent Frialit<sup>®</sup> Implant.

Dentsply-Friadent Frialit <sup>®</sup>	3.4 mm	3.8 mm	4.5 mm	5.5 mm
Ti-Base T-Series	T 800	T 805	T 810	T 820
Abutment screw	T 60	T 60	T 60	T 60

The Abutment Ti-Base H-Series is compatible with Biomet 3i Osseotite<sup>®</sup> Certain<sup>®</sup>.

Biomet 3i Osseotite <sup>®</sup> Certain <sup>®</sup>	3.4 mm	4.1 mm	5.0 mm
Ti-Base H-Series	H800	H810	H 820
Abutment screw	H 60	H 60	H 60

The Abutment Ti-Base L-Series is compatible with Straumann<sup>®</sup> Bone Level Implant.

Straumann <sup>®</sup> BoneLevel	3.3 mm	4.1 / 4.8 mm
Ti-Base L-Series	L 800	L 810
Abutment screw	L 60	L 61

The Abutment Ti-Base F-Series is compatible with Nobel Biocare NobelActive<sup>™</sup> Implant.

NobelActive <sup>™</sup>	3.5 mm	4.3 / 5.0 mm
Ti-Base F-Series	F 800	F 810
Abutment screw	F 60	F 61

Each Ti-Base is delivered with an abutment screw for fixation on the implant. The article number is the order number.

#### COMPOSITION:

Ti-Base and Abutment Screw: Ti6Al4V, medical grade 5, ASTM 136  
Scan Body: Polyether-ether-ketone, PEEK

#### CONTRAINDICATION:

The Ti-Bases of each Series can only be combined with the matching implant, e.g. the E-Series shall be combined exclusively with Replace select<sup>®</sup> implants. They cannot be combined with implants of a different implant type or manufacturer. The diameter of the Ti-base must correspond in size to the used implant in order to prevent a peri-implant tissue irritation.

The Ti-Bases are indicated for single use only. If they are used multiple times, they might damage the implants.

For fixation of the Ti-Bases on the implant, the correct torque force, recommended by the implant manufacturer, has to be considered carefully to avoid the damage of the implant-bone connection.

Mechanical treatment of the connection part of the Ti-Base will damage the correct fitting of the Ti-Base on the implant.

Nom	Abutment					
20	H-Series					
25	S-Series	T-Series	F-Series			
30	R-Series					
35	I-Series	K-Series	N-Series	E-Series	L-Series	

Mechanical treatment of the connection part of the Ti-Base will damage the correct fitting of the Ti-Base on the implant.

#### HANDLING METHOD FOR FURTHER PROCESSING:

The following instruction describes possible steps in the Laboratory for further processing to design the prosthetic components. The ceramic crown is an example for a possible aid prosthetic.

Ceramic crown: Milled with CAD/CAM-machines for zirconium oxide- or aluminum oxide- ceramic according to the anatomic form of a crown. The ceramic crown shall be grinded or polished with diamond instruments and with minimal pressure and water-cooling. The minimal thickness shall be 0.5 mm. Sharpe edges must be avoided.

#### VENEERING:

Crowns shall be veneered with appropriate ceramics before cementing onto the Ti-Base. The instructions for use of the ceramic manufacturers have to be considered.

Before cementing the Ti-Base and the ceramic crown: Sandblasting the contact surfaces with Al<sub>2</sub>O<sub>3</sub>, 50 µm, 2 bar and ensure intensive cleaning of dust and grease.

It is recommended to protect the connection part of the Ti-Base with an implant analog during handling.

#### CEMENTING:

It is recommended to cement the ceramic abutment onto the Ti-Base with Panavia<sup>®</sup> F2.0 (Kuraray) with RelayXUnicem<sup>®</sup> (3M Espe) or other equivalent cements. The instructions for use of the cements shall be followed carefully.

The Ti-Base shall be fixed onto an implant analog with the abutments screw. The head of the screw has to be covered with wax or resin. The mixed cement is applied onto the contact part of the Ti-Base. The abutment is pressed onto the Ti-Base. The final position is evaluated by slight rotation. The gap between abutment and the Ti-Base must be as small as possible. Remaining cement shall be removed immediately.

#### POLISHING:

After hardening the remaining cement shall be removed with rotating silicon instruments. The cement inside the screw channel has to be removed carefully.



**SCAN BODY INDICATIONS:**

The Scan Body can be used as auxiliary component, to determine the exact position and insertion angle of the implant. This is helpful for the further process to design the prosthetics.

**SCAN BODY:**

For the CAD/CAM scanning of the model, the Scan Body is used to indicate the position of the implant. The size of the Scan Body shall be corresponding to the original Implant system, implant diameter and Ti-Base Series. The chamfer of the Scan Body prevents the rotation of the ceramic abutment.

The Scan Body is fixed on the implant analogue with the abutment screw. After correct positioning, there is no gap visible between implant and Scan Body. Rotation of the Scan Body is impossible.

**ARTICLE NO.:**

The article number of Scan Body and Titanium Base is a combination of the code for the Series: E, I, K, N, R, S, T, H, L and F (→ X) with the code W for Scan Body and for Titanium Base.

Titanium-Base	X 800	X 805	X 810	X 820	X 830
Scan-Body	X 00W	X 05W	X 10W	X 20W	X 30W

**WARNING:**

**Safety hint:** metal dust is harmful to your health. When milling and sandblasting use a suction extraction system and a breathing mask.

**SECONDARY EFFECTS:**

Allergies to the alloy or contents of the alloy or electrochemically based reactions may very rarely occur.

**REACTIONS:**

In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

**WARRANTY:**

10 Years on the mechanical stability of the Ti-Base, if it was processed according to the Instruction for use. Whether given verbally, in writing or by practical instructions, our recommendation for use is based upon own experience and trials and can only be considered as standard values. Our products are subject to a constant further development. Therefore alternations in construction and composition are reserved.

**CLEANING, DISINFECTION AND STERILIZATION:**

The nt-trading abutments and screws of the series E, I, K, N, R, S and T are supplied in non-sterile condition. The components should be cleaned, disinfected and in specific clinical procedures and cases be sterilized, prior to use after they are received from the dental laboratory (no liability on disregard). Effective cleaning and disinfection is an indispensable requirement for effective sterilization of the abutments. Prior to sterilization, please keep implants and screws clean when handling in the laboratory and operator.

Additionally, please pay attention to legal regulations valid for your local areas as well as to the hygienic instructions of your dental practice. This applies particularly to the different guidelines regarding the inactivation of prions.

**1. Pre-disinfection** (avoidance of cross contaminations)

Place the abutments and screws in a germicidal bath\* immediately after use. Remove all residues and disassemble demountable products.

**2. Cleaning**

Please use distilled water and neutral cleaning agents\* only. The internal irrigation tube has to be cleaned with a Miller needle and must be rinsed with distilled water at the beginning and end of the exposure time using a disposable syringe (min. 10 ml). The products must be cleaned with a plastic instrument cleaning brush and then rinsed with distilled water. Please control all products after cleaning in order to avoid either damaging or corrosion. Damaged products must be replaced.

**3. Rinsing and Drying**

After removal of the products from the germicidal bath, all components must be rinsed 3 times with distilled water /e.g. Aqua purificata). Please dry all components thoroughly with a lint-free disposable cloth. For the cleaning of the internal irrigation tube oil-free compressed air is mandatory. Please re-check all parts for damage or corrosion afterwards.

**4. Disinfection**

We recommend a high level disinfectant such as, Cidex OPA (Johnson & Johnson) for disinfection of the abutments and screws.

- Soak the abutments and screws in the disinfectant solution for the required amount of time. See instructions for use of Cidex OPA.
- Remove the abutments and screws from the disinfectant solution.
- Rinse at least three times with highly purified water.
- Air dry and package the abutments and screws immediately.

**5. Sterilization:**

If no sterilization device is available in the laboratory, this information should be forwarded to the dentist so proper sterilization can occur. Please use only validated sterilization procedures for the sterilization of the abutments and screws. Other sterilization procedures must not be used.

**Reusability:**

You may only sterilize the abutments one time. In case of inadvertent contamination, you may re-sterilize one time after cleaning and disinfection.

**Steam sterilization:**

- fractionated vacuum procedure or gravity procedure (with sufficient product drying)
- steam sterilizer according to ISO 17665: 2006 or EN 13060 and EN 285 respectively or equivalent national standards
- validated according to EN ISO/ANSI AAMI 17665 (in past: EN 554/ANSI AAMI ISO 11134) (valid IQ/OQ) (commissioning and product specific performance qualification)
- sterilization time 20 minutes at 121 °C (250 °F) (listed exposure times are at sterilization temperature)

**7. Storage**

Store the sterilized parts dry and dust-free at room temperature.

\* Please observe all manufacturers guidelines for disinfection and cleaning agents with special regard to the concentration, exposure time and temperature. Only neutral disinfection solutions without chlorine, ammonia and aldehydes and with a proven effectiveness against HBV, HCV, and HIV must be used. The products have to meet the respective national regulations for disinfectants. If disinfectants containing aldehydes are used, this might lead to a possible fixation of proteins. Please use only freshly prepared solutions.

Rev. A/ 2012-03-20

**MANUFACTURER:**  
nt-trading GmbH & Co KG

Essostrasse 16  
76187 Karlsruhe  
Germany  
Tel: +49 - 721 - 91 54 71 - 60  
Fax: +49 - 721 - 91 54 71 - 61  
E-mail: info@nt-trading.com

# **Traditional 510(k) Submission**

Straumann® Variobase™ Abutments

---

## **Appendix 4 – NT-Trading Ti-Base 510(k) Summary**

510(k) Summary

K111935

FEB 17 2012

Submitter Name: NT-Trading GmbH & Co. KG  
Submitter Address: Eссоstrasse 16  
76187 Karlsruhe  
Germany

Phone Number: +49-721-915471 60  
Fax Number: +49-721-915471 61

Contact Person: Dirk Jahn

Date Prepared: June 29, 2011

Device Trade Name: Ti-Base Abutment  
2-CONnect Abutment

Common Name: Dental Abutments

Classification Name, Number & Product Code: Abutment, Implant, Dental, Endosseous  
872.3630  
NHA

Predicate Devices: (K100152) Sirona Dental Systems Sirona Dental CAD/CAM System, (K083871) Atlantis™ Straumann Bone Level Abutment, (K093483) Atlantis™ Abutment for Nobel Active Implant, (K072642) Biomet 3I Dental Abutments And Restorative Components, (K990342) synOcta® Prosthetics, (K080239) P.004 Abutments, (K072570) NobelActive™ Multi Unit Abutment

Device Description and Statement of Intended Use The Ti-Base Abutment is a premanufactured prosthetic component supplied in two parts, the abutment and screw, for fixation onto dedicated endosseous dental implants and is intended for use as an aid in prosthetic rehabilitation.

The 2-CONnect Abutment consists of 1 Abutment with screw (for fixation of abutment to the implant) and 1 titanium cap with 1 tiny screw (fixed into the hollow Abutment screw). The cap on top fits exactly to the abutment-geometry and does not have a rotation fixation, so it is easier to work with (not indicated for single crowns but strictly for bridges). The 2-CONnect is intended for use as an aid in prosthetic rehabilitation.

The NT-Trading Ti-Base and 2-CONnect is compatible with commercially available dental CAD/CAM systems, such as 3Shape, Exocad, Dental Wings. Such systems must be validated by the user.

Indication for use:

**Ti-Base Abutments:** The devices covered by this submission are abutments which are placed into a dental implant to provide support for a prosthetic restoration.

The Ti-Base abutments are intended for use to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support single and multiple tooth prosthesis, in the mandible or maxilla. The prosthesis can be cement retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.

The Ti-Base abutments are indicated for use with the following implant systems:



- Nobel Biocare® Replace Select®
- Nobel Biocare NobelActive™
- Biomet 3i® Osseotite®
- Biomet 3i® Osseotite® Certain®
- Nobel Biocare Branemark®
- Straumann® synOcta®
- Straumann® Bone Level®
- Zimmer® Tapered Screw-vent®
- Astra Tech OsseoSpeed®
- Dentsply-Friadent® Frialit®

**2-CONNECT Abutments:** 2-CONNECT abutment is indicated for use to provide support for prosthetic restorations such as bars and bridges. The 2-CONNECT abutments can be used in multiple tooth restorations. The 2-CONNECT abutment can be used together with cemented bridges and bar constructions for functional and aesthetical reconstruction.

The 2-CONNECT abutments are indicated for use with the following implant systems:

- Nobel Biocare® Replace Select®
- Straumann® synOcta®
- Straumann® BoneLevel®

Summary of Technological  
Characteristics

The proposed Ti-Base abutments and 2-CONNECT abutments are substantially equivalent to the currently cleared devices. They are substantially equivalent in intended use, material and connection interfaces to the implants are identical for each individual diameter and connection type. Comparison Demonstrating Substantial Equivalence follows at the end of this section.

**Testing Summary**

In order to demonstrate compatibility of Ti-Base and 2-CONNECT abutments to each implant system, fatigue testing was performed according to ISO 14801 Dentistry-Implants-Dynamic fatigue test for endosseous implants. Testing was performed on the abutments in this submission with the implants that they are intended to fit. See section 18.

Conclusion

The information discussed above demonstrates that the NT-Trading Ti-Base Dental Abutments and 2-CONNECT Abutments are substantially equivalent to the predicate devices.

Declarations

- This summary includes only information that is also covered in the body of the 510(k).
- This summary does not contain any puffery or unsubstantiated labeling claims.
- This summary does not contain any raw data, i.e., contains only summary data.
- This summary does not contain any trade secret or confidential commercial information.
- This summary does not contain any patient identification information.

Summary of Technical Characteristics

Feature	Ti-Base and 2-CONNECT	Sirona Dental Systems Sirona Dental CAD/CAM System	Atlantis™ Straumann Bone Level Abutment	Atlantis™ Abutment for Nobel Active Implant	Biomet 3i Dental Abutments And Restorative Components	synOcta® Prosthetics	P.004 Abutments	NobelActive™ Multi Unit Abutment
510(k) Number		K100152	K083871	K093483	K072642	K990342	K080239	K072570
Manufacturer	Ni-Trading GmbH & Co. KG	Sirona Dental Systems GmbH	Astra Tech Inc.	Astra Tech Inc.	Biomet 3i, Inc.	Straumann® USA	Straumann® Manufacturing, Inc	Nobel Biocare® AB
Classification # & Product Code	872.3630 NHA	872.3630 NHA	872.3630 NHA	872.3630 NHA	872.3630 NHA	872.3630 NHA	872.3630 NHA	872.3630 NHA
Intended Use	<p><b>Ti-Base Abutments:</b> The devices covered by this submission are abutments which are placed into a dental implant to provide support for a prosthetic restoration.</p> <p>The Ti-Base abutments are intended for use to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support single tooth prosthesis, in the mandible or maxilla. The prosthesis can be cement retained to the abutment. The abutment screw is</p>	<p>The Sirona Dental CAD/CAM System is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations. The system consists of three major parts: TiBase, InCoris mesostructure, and CAD/CAM software. Specifically, the InCoris mesostructure and TiBase components make up a two-piece abutment which is used in conjunction with endosseous dental implants to restore the function</p>	<p>The devices covered by this submission are abutments which are placed into a dental implant to provide support for a prosthetic reconstruction. The Atlantis Abutment is intended for use with an endosseous implant to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support single tooth prosthesis, in the mandible or maxilla. The prosthesis can be cement retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.</p>	<p>The Atlantis Abutment is intended for use with an endosseous implant to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support single tooth prosthesis, in the mandible or maxilla. The prosthesis can be cement retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.</p>	<p>BIOMET 3i Dental Abutments and Overdenture Bars are intended for use as an accessory to endosseous dental implants to support a prosthetic device in a partially or edentulous patient. These are intended for use to support single and multiple tooth prostheses, in the mandible or maxilla. The prostheses can be cement retained to the abutment.</p>	<p>ITI Dental implants are intended to be placed in the maxillary and/or mandibular arch to support crowns, bridges or overdentures in edentulous or partially edentulous patients. The prosthetic accessories to dental implants are used either in the process of fabricating the prosthetic restoration for the implant or as part of the prosthetic restoration.</p>	<p>Abutments are placed into dental implants to provide support for prosthetic restorations such as crowns, bridges and overdentures. Abutments can be used in single tooth replacements and multiple tooth restorations. The subject abutments are for permanent screw-retained bridges and bar-retained implant-borne dentures. Permanent copings are intended to</p>	<p>Nobel Biocare's Multi-Unit is a premanufactured prosthetic component directly connected to endosseous dental implants and is intended for use as an aid in prosthetic rehabilitation.</p>

<p>intended to secure the abutment to the endosseous implant.</p> <p>The Ti-Base abutments are indicated for use with the following implant systems:</p> <ul style="list-style-type: none"> <li>• Nobel Biocare® Replace Select®</li> <li>• Nobel Biocare™ NobelActive™</li> <li>• Biomet 3i®</li> <li>• Osseolite®</li> <li>• Biomet 3i®</li> <li>• Osseolite®</li> <li>• Certain®</li> <li>• Nobel Biocare Branemark®</li> <li>• Straumann®</li> <li>• synOcta®</li> <li>• Straumann® Bone Level®</li> <li>• Zimmer® Tapered Screw-vent®</li> <li>• Astra Tech OsseoSpeed®</li> <li>• Dentsply-Friadent® Frialit®</li> </ul> <p><b>2-CONNECT Abutments:</b> 2-CONNECT abutment is indicated for use to provide support for prosthetic restorations such as bars and bridges. The 2-CONNECT abutments can be</p>	<p>and aesthetics in the oral cavity. The InCoris mesostructure may also be used in conjunction with the Camlog Titanium base CAD/CAM (types K2244.xxxx) (K083496) in the Camlog Implant System. The CAD/CAM software is intended to design and fabricate the InCoris mesostructure. The InCoris mesostructure and TiBase two-piece abutment is compatible with the following implants systems:</p> <ul style="list-style-type: none"> <li>• Nobel Biocare Replace</li> <li>• Nobel Biocare Branemark</li> <li>• Friadent Xive</li> <li>• Biomet 3i</li> <li>• Osseolite</li> <li>• Astra Tech Osseospeed</li> <li>• Zimmer Tapered Screw-Vent</li> <li>• Straumann SynOcta</li> </ul>	<p>retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.</p>				<p>serve as a base for multi-unit bar or bridge restorations. Temporary Copings are intended to serve as a base for temporary restorations for up to 6 month. Protective Caps are intended to protect the outer configuration of the abutment and to maintain and condition the contours of the soft tissue during the healing phase for up to 6 months.</p>
---	--	---	--	--	--	--

	used in multiple tooth restorations. The 2-CONNECT abutment can be used together with cemented bridges and bar constructions for functional and aesthetical reconstruction. The 2-CONNECT abutments are indicated for use with the following implant systems:										
	<ul style="list-style-type: none"> <li>• Nobel Biocare®</li> <li>• Straumann®</li> <li>• synOcta®</li> <li>• Straumann® BoneLevel®</li> </ul>	Same	Same	Same	Same	Same	Same	Same	Same	Same	Same
Abutment Diameter min.	3.5 mm	Same	Same	Same	Same	Same	Same	Same	Same	Same	Same
Abutment Diameter max.	6.5 mm	Same	Same	Same	Same	Same	Same	Same	Same	Same	Same
Abutment Height	Ti-Base: 4 mm 2-CONNECT: 2.3 / 4.3 mm	Same	4 / 5.5 mm	Same	6.6 mm	7.0	1.5 / 6.0	1.5 / 6.0	1.5 / 6.0	1.0 / 5.5 mm	
Mode of Action	Screw-retained or cement retained	Screw-retained or cement retained	Screw-retained or cement retained	Screw-retained or cement retained	Screw-retained or cement retained	Screw-retained or cement retained	Screw-retained or cement retained	Screw-retained or cement retained	Screw-retained or cement retained	Screw retained	Screw retained
Reusable	No	No	No	No	No	No	No	No	No	No	No
Material	Ti-6Al-4V	Ti-6Al-4V	Ti-6Al-4V ELI	Ti-6Al-4V ELI	Ti-6Al-4V ELI	Ti-6Al-4V	Ti-6Al-4V	Titanium, Titanium alloy	Ti-6Al-4V	Titanium Alloy	

For the reasons stated above, we believe a determination of substantial equivalence between the Ti-Base and 2-CONNECT and these predicate devices is appropriate.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room –WO66-G609  
Silver Spring, MD 20993-0002

NT-Trading GmbH & Company AG  
C/O Mr. William Greenrose  
President  
Qserve America, Inc.  
220 River Road  
Claremont, New Hampshire 03743

FEB 17 2012

Re: K111935  
Trade/Device Name: Ti-Base for Individual milled Zirconium Abutment, 2-CONnect  
Abutment for Bridges and Bars  
Regulation Number: 21 CFR 872.3630  
Regulation Name: Endosseous Dental Implant Abutment  
Regulatory Class: II  
Product Code: NHA  
Dated: November 28, 2011  
Received: February 14, 2012

Dear Mr. Greenrose:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Anthony D. Watson, B.S., M.S., M.B.A.  
Director  
Division of Anesthesiology, General Hospital,  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K111935

Device Name: Ti-Base for individual milled Zirconium Abutment, 2-CONnect Abutment for Bridges and Bars

Indications For Use:

**Ti-Base for individual Zirconium Abutments:** The devices covered by this submission are abutments which are placed into a dental implant to provide support for a prosthetic restoration.

The Ti-Base for individual Zirconium Abutments are intended for use to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support single and multiple tooth prosthesis, in the mandible or maxilla. The prosthesis can be cement retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.

The Ti-Base abutments are indicated for use with the following implant systems:

- Nobel Biocare® Replace Select®
- Nobel Biocare NobelActive™
- Biomet 3i® Osseotite®
- Biomet 3i® Osseotite® Certain®
- Nobel Biocare Branemark®
- Straumann® synOcta®
- Straumann® Bone Level®
- Zimmer® Tapered Screw-vent®
- Astra Tech OsseoSpeed®
- Dentsply-Friadent® Frialit®

**2-CONnect Abutment for Bridges and Bars:** 2-CONnect Abutment for Bridges and Bars is indicated for use to provide support for prosthetic restorations such as bars and bridges. The 2-CONnect abutments can be used in multiple tooth restorations. The 2-CONnect abutment can be used together with cemented bridges and bar constructions for functional and aesthetical reconstruction.

The 2-CONnect abutments are indicated for use with the following implant systems:

- Nobel Biocare® Replace Select®
- Straumann® synOcta®
- Straumann® BoneLevel®

Prescription Use  (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
(Division Sign-Off)

Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K111935

# **Traditional 510(k) Submission**

Straumann® Variobase™ Abutments

---

## **Appendix 5 – Biocompatibility Test Report for Straumann® Variobase™ Abutments**











































# COVER SHEET MEMORANDUM

Food and Drug Administration  
Office of Device Evaluation &  
Office of In Vitro Diagnostics and  
Radiological Health

**NOTE: This form is REQUIRED for holds and for final decisions.**

Reviewer Name Michael Mendelson

510(k) Number K132219

**Please list CTS decision code:** SE - Substantially Equivalent

Hold (Additional Information or Telephone Hold) Hold Date

Final Decision (SE, SE with Limitations, NSE (select code below), Withdrawn, etc.)

Incomplete Response - Convert Supplement to Amendment (attach email sent to firm)

Add to File

(review staff should follow the instructions and complete the memo/routing sheet at:  
[http://eroom.fda.gov/eRoom/CDRH3/CDRHPreMarketNotification510kProgram/0\\_3bba7](http://eroom.fda.gov/eRoom/CDRH3/CDRHPreMarketNotification510kProgram/0_3bba7). DCC should refer to that documentation for the close-out code and mail any provided letter.)

The remainder of this form must be filled out for close-outs only

<b>Class:</b>	II
<b>Regulation Number:</b>	872.3630
<b>Product Code:</b>	NHA
<b>Additional Product Codes:</b>	



Please complete the following for a final clearance decision (i.e, SE, SE with Limitations, etc.)	YES	NO
Indications for Use Page (Attach IFU)	X	
510(k) Summary or 510(k) Statement (Attach Summary)	X	
Truthful and Accurate Statement (Must be present for a Final Decision)	X	
Is the device Class III?		X
Is this a combination product?		X
Is this device intended for pediatric use only?		X
Is this a prescription device? (If both prescription & OTC, check both boxes.)	X	
Is clinical data necessary to support the review of this 510(k)?		X
For United States based clinical studies only, did the application include a completed Form FDA 3674, Certification with Requirements of ClinicalTrials.gov Data Bank? (If study was conducted in the United States and Form FDA 3674 was not included or was incomplete, then applicant must be contacted to obtain completed form.)		
Does this device include an Animal Tissue Source?		X
All Pediatric Patients age <= 21		X

DCC-  
2/21

Neonate/Newborn (Birth to 28 days)		×
Infant (29 days to < 2 years)		×
Child (2 years to <12 years)		×
Adolescent (12 years to <18 years)		×
Transitional Adolescent A (18 years to <21 years); Special considerations are being given to this group, different from adults age >= 21 (different device design or testing, different protocol procedures, etc.)		×
Transitional Adolescent B (18 years to <21 years); No special considerations compared to adults >= 21 years)		×
Is this device subject to the Tracking Regulation? ( <a href="#">Medical Device Tracking Guidance</a> )		×

**Digital Signature Concurrence Table**

(Not all signatures may be required)

Branch Chief Sign-Off	 <p>Mary S. Runner -S  2014.02.21  10:44:07 -05'00'</p>
Division Sign-Off	 <p>Kwame O. Ulmer -S  2014.02.21 14:01:22 -05'00'</p>



K132219/S 001

FDA CDRH DMC

AUG 07 2013

Received

August 6, 2013

U.S. Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center – W066-G609  
10903 New Hampshire Avenue  
Silver Spring, Maryland 20993-0002

**Subject: K132219 Response to RTA Letter dated August 5, 2013**

Dear Sir or Madam:

On behalf of Institut Straumann AG, Straumann USA, LLC submits this response to the Refuse to Accept Letter received on 05-Aug-2013 for K132219 Traditional 510(k) Premarket Notification for the Straumann® Variobase™ Abutments in accordance with 21 CFR 807.81.

(b)(4) Trade Secret Process - Product Specs



(b)(4) Trade Secret Process - Product Specs



(b)(4) Trade Secret Process - Product Specs



(b)(4) Trade Secret Process - Product Specs



(b)(4) Trade Secret Process - Product Specs


A large, solid black rectangular redaction box covers the majority of the page's content, starting below the header and extending nearly to the bottom. The text '(b)(4) Trade Secret Process - Product Specs' is printed in red at the top left corner of this redacted area.



(b)(4) Trade Secret Process - Product Specs



Sincerely,



Jennifer M. Jackson, MS  
Senior Regulatory Affairs Project Manager

Enclosures



August 6, 2013

U.S. Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center – W066-G609  
10903 New Hampshire Avenue  
Silver Spring, Maryland 20993-0002

**Subject: K132219 Response to RTA Letter dated August 5, 2013**

Dear Sir or Madam:

On behalf of Institut Straumann AG, Straumann USA, LLC submits this response to the Refuse to Accept Letter received on 05-Aug-2013 for K132219 Traditional 510(k) Premarket Notification for the Straumann® Variobase™ Abutments in accordance with 21 CFR 807.81.

(b)(4) Trade Secret Process - Product Specs



(b)(4) Trade Secret Process - Product Specs

A large, solid black rectangular redaction box covers the majority of the page content, starting below the header and ending above the footer. The text '(b)(4) Trade Secret Process - Product Specs' is written in red at the top left corner of this redacted area.

(b)(4) Trade Secret Process - Product Specs

A large, solid black rectangular redaction box covers the majority of the page content, starting below the header and ending above the footer. The text '(b)(4) Trade Secret Process - Product Specs' is printed in red at the top left corner of this redacted area.

(b)(4) Trade Secret Process - Product Specs

A large, solid black rectangular redaction box covers the majority of the page content, starting below the header and ending above the footer. The text '(b)(4) Trade Secret Process - Product Specs' is printed in red at the top left corner of this redacted area.

(b)(4) Trade Secret Process - Product Specs

A large, solid black rectangular redaction box covers the majority of the page content, starting below the header and ending above the footer. The text '(b)(4) Trade Secret Process - Product Specs' is printed in red at the top left corner of this redacted area.



Two copies of the response are enclosed. The second copy is being provided in PDF format on CD and is an exact duplicate of the hardcopy with the exception that the eCopy does not contain original signatures.

We trust that the foregoing information will be sufficient to permit FDA to make a finding of substantial equivalence for the proposed Straumann Dental Implant System to the currently marketed devices as presented in this premarket notification.

Please address any questions regarding this 510(k) Premarket Notification to the undersigned.

Sincerely,

A handwritten signature in blue ink that reads "Jennifer M. Jackson".

Jennifer M. Jackson, MS  
Senior Regulatory Affairs Project Manager

Enclosures

K132219/S2

September 13, 2013

U.S. Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center – W066-G609  
10903 New Hampshire Avenue  
Silver Spring, Maryland 20993-0002

FDA CDRH DMC

SEP 16 2013

Received

**Subject: K132219 Response to RTA Letter dated August 16, 2013**

Dear Sir or Madam:

On behalf of Institut Straumann AG, Straumann USA, LLC submits this response to the Refuse to Accept Letter received on 16-Aug-2013 for K132219 Traditional 510(k) Premarket Notification for the Straumann® Variobase™ Abutments in accordance with 21 CFR 807.81.

(b)(4) Trade Secret Process - Product Specs



41




(b)(4) Trade Secret Process - Product Specs

A large, solid black rectangular redaction box covers the majority of the page's content, starting below the header and extending nearly to the bottom. The text '(b)(4) Trade Secret Process - Product Specs' is printed in red at the top left corner of this redacted area.

(b)(4) Trade Secret Process - Product Specs

A large, solid black rectangular redaction box covers the majority of the page's content, starting below the redaction code and ending above the signature.

Sincerely,

A handwritten signature in blue ink that reads 'Jennifer M. Jackson'.

Jennifer M. Jackson, MS  
Senior Regulatory Affairs Project Manager

Enclosures



September 13, 2013

U.S. Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center – W066-G609  
10903 New Hampshire Avenue  
Silver Spring, Maryland 20993-0002

**Subject: K132219 Response to RTA Letter dated August 16, 2013**

Dear Sir or Madam:

On behalf of Institut Straumann AG, Straumann USA, LLC submits this response to the Refuse to Accept Letter received on 16-Aug-2013 for K132219 Traditional 510(k) Premarket Notification for the Straumann® Variobase™ Abutments in accordance with 21 CFR 807.81.

(b)(4) Trade Secret Process - Product Specs



(b)(4) Trade Secret Process - Product Specs

A large, solid black rectangular redaction box covers the majority of the page content, starting below the header and ending above the footer.

(b)(4) TS/CCI

Information in this premarket notification is considered proprietary or trade secret or confidential commercial information. The company requests that all such information not be disclosed pursuant to 18 U.S.C. §1905, 5 U.S.C. §552, 21 U.S.C. §331(j), and all other applicable laws and regulations.

We trust that the foregoing information will be sufficient to permit FDA to make a finding of substantial equivalence for the proposed Straumann® Variobase™ Abutments to the currently marketed devices as presented in this premarket notification.

Please address any questions regarding this 510(k) Premarket Notification to the undersigned.

Sincerely,



Jennifer M. Jackson, MS  
Senior Regulatory Affairs Project Manager

Enclosures

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Table of Contents

---

**Table of Contents**

1	Medical Device User Fee Cover Sheet (Form FDA 3601) .....	1-1
2	CDRH Premarket Review Submission Cover Sheet.....	2-1
3	510(k) Cover Letter .....	3-1
4	Indications for Use Statement.....	4-1
5	510(k) Summary .....	5-1
5.1	Submitter's Contact Information .....	5-1
5.2	Name of the Device .....	5-1
5.3	Predicate Device(s) .....	5-1
5.4	Device Description.....	5-1
5.5	Intended Use .....	5-1
5.6	Technological Characteristics.....	5-2
5.7	Performance Testing .....	5-2
5.8	Conclusion.....	5-2
6	Truthful and Accuracy Statement.....	6-1
7	Class III Summary and Certification.....	7-1
8	Financial Certification or Disclosure Statement .....	8-1
9	Declarations of Conformity and Summary Reports.....	9-1


(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Table of Contents

(b)(4) Trade Secret Process - Product Specs



10 Device Description .....10-1

(b)(4) Trade Secret Process - Product Specs



11 Substantial Equivalence Discussion.....11-1


12 Proposed Labeling .....12-1

    12.1 Package Label .....12-1

    12.2 Proposed Package Insert/Instructions for Use for Straumann® Variobase™  
    Abutments .....12-2

13 Sterilization and Shelf Life .....13-1

(b)(4) Trade Secret Process - Product Specs



14 Biocompatibility .....14-1

(b) [Redacted] .....14-1

(4)

15 Software .....15-1

16 Electromagnetic Compatibility and Electrical Safety .....16-1

17 Performance Testing – Bench .....17-1





**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Table of Contents

(	
b	
)	
(	
4	
)	
T	
r	
a	
d	
e	
S	1

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Table of Contents

---

**Table of Tables**

(b)(4) Trade Secret Process - Product Specs

[Redacted Table of Contents]

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Table of Contents

---

**Table of Figures**

[REDACTED]	(
[REDACTED]	b
[REDACTED]	)
[REDACTED]	4
[REDACTED]	)
[REDACTED]	T
[REDACTED]	r
[REDACTED]	a
[REDACTED]	d
[REDACTED]	e

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

Medical Device User Fee Cover Sheet (Form FDA 3601)

---

### 1 Medical Device User Fee Cover Sheet (Form FDA 3601)

Payment Identification Number: (b)(4) Trade Secret  
Process Product

The Medical Device User Fee Cover Sheet (Form FDA 3601) begins on the next page.

Form Approved OMB No. 0910-0511 Expiration Date April 30, 2016. See Instructions for OMB Statement.

DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION <b>MEDICAL DEVICE USER FEE COVER SHEET</b>		PAYMENT IDENTIFICATION NUMBER: <b>(b)(4) Trade</b> Write the Payment Identification number on your check.
A completed cover sheet must accompany each original application or supplement subject to fees. If payment is sent by U.S. mail or courier, please include a copy of this completed form with payment. Payment and mailing instructions can be found at: <a href="http://www.fda.gov/oc/mdufma/coversheet.html">http://www.fda.gov/oc/mdufma/coversheet.html</a>		
1. COMPANY NAME AND ADDRESS (include name, street address, city state, country, and post office code)  STRAUMANN USA 60 MINUTEMAN ROAD ANDOVER MA 01810 US  1.1 EMPLOYER IDENTIFICATION NUMBER (EIN) <b>(b)(4)</b>	2. CONTACT NAME Jennifer Jackson 2.1 E-MAIL ADDRESS jennifer.jackson@straumann.com 2.2 TELEPHONE NUMBER (include Area code) 978-747-2509 2.3 FACSIMILE (FAX) NUMBER (Include Area code) 978-747-0023	
3. TYPE OF PREMARKET APPLICATION (Select one of the following in each column; if you are unsure, please refer to the application descriptions at the following web site: <a href="http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm345263.htm">http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm345263.htm</a> ) <u>Select an application type:</u> <input checked="" type="checkbox"/> Premarket notification(510(k)); except for third party <input type="checkbox"/> 513(g) Request for Information <input type="checkbox"/> Biologics License Application (BLA) <input type="checkbox"/> Premarket Approval Application (PMA) <input type="checkbox"/> Modular PMA <input type="checkbox"/> Product Development Protocol (PDP) <input type="checkbox"/> Premarket Report (PMR) <input type="checkbox"/> 30-Day Notice  3.1 Select a center <input checked="" type="checkbox"/> CDRH <input type="checkbox"/> CBER <u>3.2 Select one of the types below</u> <input checked="" type="checkbox"/> Original Application <u>Supplement Types:</u> <input type="checkbox"/> Efficacy (BLA) <input type="checkbox"/> Panel Track (PMA, PMR, PDP) <input type="checkbox"/> Real-Time (PMA, PMR, PDP) <input type="checkbox"/> 180-day (PMA, PMR, PDP)		
4. ARE YOU A SMALL BUSINESS? (See the instructions for more information on determining this status) <input type="checkbox"/> YES, I meet the small business criteria and have submitted the required qualifying documents to FDA <input checked="" type="checkbox"/> NO, I am not a small business 4.1 If Yes, please enter your Small Business Decision Number:		
5. FDA WILL NOT ACCEPT YOUR SUBMISSION IF YOUR COMPANY HAS NOT PAID AN ESTABLISHMENT REGISTRATION FEE THAT IS DUE TO FDA. HAS YOUR COMPANY PAID ALL ESTABLISHMENT REGISTRATION FEES THAT ARE DUE TO FDA? <input checked="" type="checkbox"/> YES (All of our establishments have registered and paid the fee, or this is our first device, and we will register and pay the fee within 30 days of FDA's approval/clearance of this device.) <input type="checkbox"/> NO (If "NO," FDA will not accept your submission until you have paid all fees due to FDA. This submission will not be processed; see <a href="http://www.fda.gov/cdrh/mdufma">http://www.fda.gov/cdrh/mdufma</a> for additional information)		
6. IS THIS PREMARKET APPLICATION COVERED BY ANY OF THE FOLLOWING USER FEE EXCEPTIONS? IF SO, CHECK THE APPLICABLE EXCEPTION. <input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates <input type="checkbox"/> This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only <input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population <input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially		
7. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION FOR WHICH FEES WERE WAIVED DUE TO SOLE USE IN A PEDIATRIC POPULATION THAT NOW PROPOSES CONDITION OF USE FOR ANY ADULT POPULATION? (If so, the application is subject to the fee that applies for an original premarket approval application (PMA)). <input type="checkbox"/> YES <input checked="" type="checkbox"/> NO		
PAPERWORK REDUCTION ACT STATEMENT Public reporting burden for this collection of information is estimated to average 18 minutes per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to the address below.  Department of Health and Human Services, Food and Drug Administration, Office of Chief Information Officer, 1350 Piccard Drive, 4th Floor Rockville, MD 20850 [Please do NOT return this form to the above address, except as it pertains to comments on the burden estimate.]		
8. USER FEE PAYMENT AMOUNT SUBMITTED FOR THIS PREMARKET APPLICATION <b>(b)(4)</b>		03-Jul-2013

Form FDA 3601 (01/2007)

["Close Window"](#) [Print Cover sheet](#)

# **Traditional 510(k) Submission**

**Straumann<sup>®</sup> Variobase<sup>™</sup> Abutments**

CDRH Premarket Review Submission Cover Sheet

---

## **2 CDRH Premarket Review Submission Cover Sheet**

The CDRH Premarket Review Submission Cover Sheet begins on the next page.

**CDRH PREMARKET REVIEW SUBMISSION COVER SHEET**

Date of Submission 09-13-2013	User Fee Payment ID Number <b>(b)(4) Trade</b>	FDA Submission Document Number (if known) K132219
----------------------------------	---	--

**SECTION A TYPE OF SUBMISSION**

<b>PMA</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Premarket Report <input type="checkbox"/> Modular Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment <input type="checkbox"/> Licensing Agreement	<b>PMA &amp; HDE Supplement</b> <input type="checkbox"/> Regular (180 day) <input type="checkbox"/> Special <input type="checkbox"/> Panel Track (PMA Only) <input type="checkbox"/> 30-day Supplement <input type="checkbox"/> 30-day Notice <input type="checkbox"/> 135-day Supplement <input type="checkbox"/> Real-time Review <input type="checkbox"/> Amendment to PMA & HDE Supplement <input type="checkbox"/> Other	<b>PDP</b> <input type="checkbox"/> Original PDP <input type="checkbox"/> Notice of Completion <input type="checkbox"/> Amendment to PDP	<b>510(k)</b> <input checked="" type="checkbox"/> Original Submission: <input checked="" type="checkbox"/> Traditional <input type="checkbox"/> Special <input type="checkbox"/> Abbreviated (Complete section I, Page 5) <input type="checkbox"/> Additional Information <input type="checkbox"/> Third Party	<b>Request for Feedback</b> <input type="checkbox"/> Pre-Submission <input type="checkbox"/> Informational Meeting <input type="checkbox"/> Submission Issue Meeting <input type="checkbox"/> Day 100 Meeting <input type="checkbox"/> Agreement Meeting <input type="checkbox"/> Determination Meeting <input type="checkbox"/> Study Risk Determination <input type="checkbox"/> Other (specify):
<b>IDE</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement	<b>Humanitarian Device Exemption (HDE)</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment	<b>Class II Exemption Petition</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information	<b>Evaluation of Automatic Class III Designation (De Novo)</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information	<b>Other Submission</b> <input type="checkbox"/> 513(g) <input type="checkbox"/> Other (describe submission):

Have you used or cited Standards in your submission?  Yes  No (If Yes, please complete Section I, Page 5)

**SECTION B SUBMITTER, APPLICANT OR SPONSOR**

Company / Institution Name Straumann USA, LLC	Establishment Registration Number (if known) 1222315		
Division Name (if applicable)	Phone Number (including area code) 978-747-2509		
Street Address 60 Minuteman Road	FAX Number (including area code) 978-747-0023		
City Andover	State / Province MA	ZIP/Postal Code 01810	Country USA
Contact Name Jennifer M. Jackson, MS			
Contact Title Senior Regulatory Affairs Project Manager		Contact E-mail Address jennifer.jackson@straumann.com	

**SECTION C APPLICATION CORRESPONDENT (e.g., consultant, if different from above)**

Company / Institution Name			
Division Name (if applicable)		Phone Number (including area code)	
Street Address		FAX Number (including area code)	
City	State / Province	ZIP Code	Country
Contact Name			
Contact Title		Contact E-mail Address	

**SECTION D1 REASON FOR APPLICATION - PMA, PDP, OR HDE**

<input type="checkbox"/> New Device <input type="checkbox"/> Withdrawal <input type="checkbox"/> Additional or Expanded Indications <input type="checkbox"/> Request for Extension <input type="checkbox"/> Post-approval Study Protocol <input type="checkbox"/> Request for Applicant Hold <input type="checkbox"/> Request for Removal of Applicant Hold <input type="checkbox"/> Request to Remove or Add Manufacturing Site	<input type="checkbox"/> Change in design, component, or specification: <input type="checkbox"/> Software / Hardware <input type="checkbox"/> Color Additive <input type="checkbox"/> Material <input type="checkbox"/> Specifications <input type="checkbox"/> Other ( <i>specify below</i> )	<input type="checkbox"/> Location change: <input type="checkbox"/> Manufacturer <input type="checkbox"/> Sterilizer <input type="checkbox"/> Packager
<input type="checkbox"/> Process change: <input type="checkbox"/> Manufacturing <input type="checkbox"/> Packaging <input type="checkbox"/> Sterilization <input type="checkbox"/> Other ( <i>specify below</i> )	<input type="checkbox"/> Labeling change: <input type="checkbox"/> Indications <input type="checkbox"/> Instructions <input type="checkbox"/> Performance Characteristics <input type="checkbox"/> Shelf Life <input type="checkbox"/> Trade Name <input type="checkbox"/> Other ( <i>specify below</i> )	<input type="checkbox"/> Report Submission: <input type="checkbox"/> Annual or Periodic <input type="checkbox"/> Post-approval Study <input type="checkbox"/> Adverse Reaction <input type="checkbox"/> Device Defect <input type="checkbox"/> Amendment
<input type="checkbox"/> Response to FDA correspondence:		<input type="checkbox"/> Change in Ownership <input type="checkbox"/> Change in Correspondent <input type="checkbox"/> Change of Applicant Address

Other Reason (*specify*):

**SECTION D2 REASON FOR APPLICATION - IDE**

<input type="checkbox"/> New Device <input type="checkbox"/> New Indication <input type="checkbox"/> Addition of Institution <input type="checkbox"/> Expansion / Extension of Study <input type="checkbox"/> IRB Certification <input type="checkbox"/> Termination of Study <input type="checkbox"/> Withdrawal of Application <input type="checkbox"/> Unanticipated Adverse Effect <input type="checkbox"/> Notification of Emergency Use <input type="checkbox"/> Compassionate Use Request <input type="checkbox"/> Treatment IDE <input type="checkbox"/> Continued Access	<input type="checkbox"/> Change in: <input type="checkbox"/> Correspondent / Applicant <input type="checkbox"/> Design / Device <input type="checkbox"/> Informed Consent <input type="checkbox"/> Manufacturer <input type="checkbox"/> Manufacturing Process <input type="checkbox"/> Protocol - Feasibility <input type="checkbox"/> Protocol - Other <input type="checkbox"/> Sponsor	<input type="checkbox"/> Response to FDA Letter Concerning: <input type="checkbox"/> Conditional Approval <input type="checkbox"/> Deemed Approved <input type="checkbox"/> Deficient Final Report <input type="checkbox"/> Deficient Progress Report <input type="checkbox"/> Deficient Investigator Report <input type="checkbox"/> Disapproval <input type="checkbox"/> Request Extension of Time to Respond to FDA <input type="checkbox"/> Request Meeting <input type="checkbox"/> Request Hearing
<input type="checkbox"/> Report submission: <input type="checkbox"/> Current Investigator <input type="checkbox"/> Annual Progress Report <input type="checkbox"/> Site Waiver Report <input type="checkbox"/> Final		

Other Reason (*specify*):

**SECTION D3 REASON FOR SUBMISSION - 510(k)**

<input checked="" type="checkbox"/> New Device	<input type="checkbox"/> Additional or Expanded Indications	<input type="checkbox"/> Change in Technology
--	---	---

Other Reason (*specify*):



**SECTION E ADDITIONAL INFORMATION ON 510(K) SUBMISSIONS**

Product codes of devices to which substantial equivalence is claimed								Summary of, or statement concerning, safety and effectiveness information	
1	NHA	2		3		4		<input checked="" type="checkbox"/> 510 (k) summary attached <input type="checkbox"/> 510 (k) statement	
5		6		7		8			

Information on devices to which substantial equivalence is claimed (if known)

	510(k) Number	Trade or Proprietary or Model Name	Manufacturer
1	K120822	Straumann CARES Variobase Abutment NNC, RN, WN, NC, RC	Straumann
2	K111935	Ti-Base Abutment	NT-Trading GmbH & Co. KG
3			
4			
5			
6			

**SECTION F PRODUCT INFORMATION - APPLICATION TO ALL APPLICATIONS**

Common or usual name or classification name

Endosseous dental implant abutment

(b)(4) Trade Secret Process - Product Specs

FDA document numbers of all prior related submissions (regardless of outcome)

1	None	2		3		4		5		6	
7		8		9		10		11		12	

Data Included in Submission

- Laboratory Testing
  Animal Trials
  Human Trials

**SECTION G PRODUCT CLASSIFICATION - APPLICATION TO ALL APPLICATIONS**

Product Code	C.F.R. Section (if applicable)	Device Class
NHA	21 CFR 872.3630	<input type="checkbox"/> Class I <input checked="" type="checkbox"/> Class II <input type="checkbox"/> Class III <input type="checkbox"/> Unclassified
Classification Panel		
Dental		

Indications (from labeling)

The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.

**Note:** Submission of the information entered in Section H does not affect the need to submit device establishment registration.

FDA Document Number (if known)

**SECTION H MANUFACTURING / PACKAGING / STERILIZATION SITES RELATING TO A SUBMISSION**

<input checked="" type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete		Facility Establishment Identifier (FEI) Number		<input checked="" type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Contract Manufacturer <input type="checkbox"/> Repackager / Relabeler	
Company / Institution Name Institut Straumann AG			Establishment Registration Number 9613348		
Division Name (if applicable)			Phone Number (including area code) 978-747-2509		
Street Address Peter Merian-Weg 12			FAX Number (including area code) 978-747-0023		
City Basel		State / Province	ZIP Code CH-4052	Country Switzerland	
Contact Name Jennifer M. Jackson, MS		Contact Title Senior Regulatory Affairs Project Manager		Contact E-mail Address jennifer.jackson@straumann.com	

<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete		Facility Establishment Identifier (FEI) Number		<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Contract Manufacturer <input type="checkbox"/> Repackager / Relabeler	
Company / Institution Name			Establishment Registration Number		
Division Name (if applicable)			Phone Number (including area code)		
Street Address			FAX Number (including area code)		
City		State / Province	ZIP Code	Country	
Contact Name		Contact Title		Contact E-mail Address	

<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete		Facility Establishment Identifier (FEI) Number		<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Contract Manufacturer <input type="checkbox"/> Repackager / Relabeler	
Company / Institution Name			Establishment Registration Number		
Division Name (if applicable)			Phone Number (including area code)		
Street Address			FAX Number (including area code)		
City		State / Province	ZIP Code	Country	
Contact Name		Contact Title		Contact E-mail Address	

## SECTION I

## UTILIZATION OF STANDARDS

**Note:** Complete this section if your application or submission cites standards or includes a "Declaration of Conformity to a Recognized Standard" statement.

	Standards No.	Standards Organization	Standards Title	Version	Date
1	ISO 10993-1:2009	ISO	Biological evaluation of medical devices -- Part 1: Evaluation and testing within a risk management process	2009	10/04/2010
2	ISO 10993-5:2009	ISO	Biological evaluation of medical devices -- Part 5: Tests for in vitro cytotoxicity	2009	05/05/2010
3	ISO 10993-12:2007	ISO	Biological evaluation of medical devices -- Part 12: Sample preparation and reference materials	2007	10/04/2010
4	ISO 10993-18:2005	ISO	Biological evaluation of medical devices -- Part 18: Chemical characterization of materials	2005	07/01/2005
5	ISO 5832-11:1994	ISO	Implants for surgery - Metallic materials -- Part 11: Wrought titanium 6-aluminum 7-niobium alloy	1994	09/09/2008
6	ISO 17665-1:2006	ISO	Sterilization of health care products -- Moist heat -- Part 1: Requirements for the development, validation, and routine control of a sterilization process for medical devices	2006	09/08/2009
7	ISO/TS 17665-2:2009	ISO	Sterilization of health care products -- Moist heat -- Part 2: Guidance on the application of ISO 17665-1	2009	01/05/2009

**Please include any additional standards to be cited on a separate page.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 0.5 hour per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
1350 Piccard Drive, Room 400  
Rockville, MD 20850

*An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number.*

**SECTION I**

**UTILIZATION OF STANDARDS**

**Note:** Complete this section if your application or submission cites standards or includes a "Declaration of Conformity to a Recognized Standard" statement.

	Standards No.	Standards Organization	Standards Title	Version	Date
1	ISO 17664:2004	ISO	Sterilization of medical devices -- Information to be provided by the manufacturer for the processing of resterilizable medical devices	2004	10-12-2012
2	ISO 14801:2007	ISO	Dentistry -- Implants -- Dynamic fatigue test for endosseous dental implants	2007	04-25-2012
3					
4					
5					
6					
7					

**Please include any additional standards to be cited on a separate page.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 0.5 hour per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
 Food and Drug Administration  
 Office of Chief Information Officer  
 Paperwork Reduction Act (PRA) Staff  
 1350 Piccard Drive, Room 400  
 Rockville, MD 20850

*An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number.*

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
510(k) Cover Letter

---

**3 510(k) Cover Letter**

The 510(k) Cover Letter begins on the next page.



September 13, 2013

U.S. Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center – W066-G609  
10903 New Hampshire Avenue  
Silver Spring, Maryland 20993-0002

**Subject: Traditional 510(k) Premarket Notification: Straumann® Variobase™ Abutments**

Dear Sir or Madam:

On behalf of Institut Straumann AG, Straumann USA, LLC submits this Traditional 510(k) Premarket Notification for the Straumann® Variobase™ Abutments in accordance with 21 CFR 807.81. The intended use and fundamental operating principles of the proposed devices are substantially equivalent to previously cleared devices as detailed in this premarket submission.

**Submitter:**

Straumann USA, LLC (on behalf of Institut Straumann AG)  
60 Minuteman Road  
Andover, MA 01810

**Primary Contact:**

Jennifer M. Jackson, MS  
Senior Regulatory Affairs Project Manager  
Straumann USA, LLC  
60 Minuteman Road  
Andover, MA 01810  
Telephone: 800-448-8168 x2509  
Fax: 978-747-0023

**Classification Name of Device:**

Classification Name:	Abutment, Implant, Dental, Endosseous
Device Product Code:	NHA
Product Classification:	Class II
Panel:	Dental
Regulation Number:	§872.3630
Prior Related Submissions:	None



**Design and Use of the Device:**

Question	YES	NO
Is the device intended for prescription use (21 CFR 801 Subpart D)?	X	
Is the device intended for over-the counter use (21 CFR 807 Subpart C)?		X
Does the device contain components derived from a tissue or other biologic source?		X
Is the device provided sterile?		X
Is the device intended for single use?	X	
Is the device a reprocessed single use device?		X
If yes, does this device type require reprocessed validation data?		
Does the device contain a drug?		X
Does the device contain a biologic?		X
Does the device use software?		X
Does the submission include clinical information?		X
Is the device implanted?		X

The term "substantially equivalent" as used herein is intended to be a determination of substantial equivalence from an FDA-regulatory point of view under the Federal Food, Drug and Cosmetic Act and relates to the fact that the product can be marketed without premarket approval or reclassification. These products may be considered distinct from a patent point of view. The term "substantially equivalent" is not applicable to, and does not diminish, any patent claim related to this product or the technology used to manufacture the product.

Information in this premarket notification is considered proprietary or trade secret or confidential commercial information. The company requests that all such information not be disclosed pursuant to 18 U.S.C. §1905, 5 U.S.C. §552, 21 U.S.C. §331(j), and all other applicable laws and regulations.

Two copies of the Traditional 510(k) Premarket Notification are enclosed. The second copy is being provided in PDF format on CD and is an exact duplicate of the hardcopy with the exception that the eCopy does not contain original signatures. Further, in accordance with the Medical Device User Fee and Modernization Act of 2001 ("MDUFMA"), Straumann USA, LLC has submitted the appropriate application fees. A copy of the User Fee Cover Sheet is provided with the enclosed 510(k) Premarket Notification.

We trust that the foregoing information will be sufficient to permit FDA to make a finding of substantial equivalence for the proposed Straumann® Variobase™ Abutments to the currently marketed devices as presented in this premarket notification.

Please address any questions regarding this 510(k) Premarket Notification to the undersigned.

Sincerely,

Jennifer M. Jackson, MS  
Senior Regulatory Affairs Project Manager

Enclosures

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Indications for Use Statement

---

**4 Indications for Use Statement**

The Indications for Use Statement associated with this 510(k) is located on the following page in the required format.



## Indications for Use

510(k) Number (if known):

Device Name: Straumann® Variobase™ Abutments

Indications for Use:

The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF  
NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
510(k) Summary

---

## **5 510(k) Summary**

### **5.1 Submitter's Contact Information**

Straumann USA, LLC (on behalf of Institut Straumann AG)

60 Minuteman Road

Andover, MA 01810

Phone Number: 1-978-747-2509

Fax Number: 1-978-747-0023

Contact Person: Jennifer M. Jackson, MS

Date of Submission: 13-Sep-2013

### **5.2 Name of the Device**

Trade Name: Straumann® Variobase™ Abutments

Common Name: Dental Implant Abutment

Classification Name: Abutment, Implant, Dental, Endosseous

Regulation Number: §872.3630

### **5.3 Predicate Device(s)**

K120822 – Straumann CARES Variobase Abutment NNC, RN, WN, NC, RC

K111935 – Ti-Base Abutment (NT-Trading GmbH & Co. KG)

### **5.4 Device Description**

The Straumann® Variobase™ Abutments are pre-manufactured (stock) abutments, sometimes referred to as “bonding bases”. Straumann® Variobase™ Abutments are available to fit Straumann® dental implant platforms NNC (Narrow Neck CrossFit®), RN (Regular Neck), WN (Wide Neck), NC (Narrow CrossFit®), and RC (Regular CrossFit®).

### **5.5 Intended Use**

Straumann® Variobase™ Abutments are pre-manufactured (stock) abutments intended to be placed onto Straumann dental implants to provide support for customized prosthetic

# **Traditional 510(k) Submission**

## **Straumann® Variobase™ Abutments**

### 510(k) Summary

---

restorations such as crowns and bridges. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.

#### **5.6 Technological Characteristics**

Straumann® Variobase™ Abutments are pre-manufactured (stock) abutments made from a titanium-aluminum-niobium alloy.

#### **5.7 Performance Testing**

The material used in the manufacture of Straumann® Variobase™ Abutments is a titanium-aluminum-niobium alloy which meets the requirements of ISO 5832-11. Bench testing was performed to evaluate the fatigue load limits of the proposed Straumann® Variobase™ Abutments. Dynamic fatigue tests were conducted in accordance with the FDA guidance document *“Guidance for Industry and FDA Staff – Class II Special Controls Guidance Document: Root-form Endosseous Dental Implants and Endosseous Dental Abutments”*.

#### **5.8 Conclusion**

The documentation submitted in this premarket notification demonstrates that the Straumann® Variobase™ Abutments are substantially equivalent to the predicate devices.


**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Truthful and Accuracy Statement

---

**6 Truthful and Accuracy Statement**

**As Required by 21 CFR 807.87(k)**

I certify that, in my capacity as Senior Regulatory Affairs Project Manager of Straumann USA, LLC, I believe to the best of my knowledge, that all data and information submitted in the premarket notification are truthful and accurate and that no material fact has been omitted.

  
(Signature)

Jennifer M. Jackson, MS  
(Typed Name)

13-Sep-2013  
(Date)

K132219  
(Premarket Notification Number)

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Class III Summary and Certification

---

**7 Class III Summary and Certification**

This section is not applicable as the subject devices have been determined to be Class II per 21 CFR 872.3630.

## **Traditional 510(k) Submission**

**Straumann® Variobase™ Abutments**

Financial Certification or Disclosure Statement

---

### **8 Financial Certification or Disclosure Statement**

This section is not applicable as there is no clinical data being submitted to support this premarket notification.

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

Declarations of Conformity and Summary Reports

---

### 9 Declarations of Conformity and Summary Reports

(b)(4) Trade Secret Process - Product Specs



# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Declarations of Conformity and Summary Reports

(b)(4) Trade Secret Process - Product Specs





# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Declarations of Conformity and Summary Reports

(b)(4) Trade Secret Process Product Specs		(b)(4) Trade Secret Process - Product Specs
A.	Applicable recognized consensus standard:	
B.	Requirements met?	Yes
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change
D.	Requirements that were not applicable to the device:	N/A
E.	Deviations from each applicable standard that were met including justification:	N/A
F.	Differences exist, if any, between the tested device to be marketed:	No
G.	Test laboratory:	(b)(4) Trade Secret Process - Product Specs

Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION  
 Traditional                       Special                       Abbreviated

STANDARD TITLE<sup>1</sup>  
**(b)(4) Trade Secret Process - Product Specs**

Please answer the following questions	Yes	No
Is this standard recognized by FDA <sup>2</sup> ? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
FDA Recognition number <sup>3</sup> .....	<b>#2-156</b>	
Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Is a summary report <sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
If no, complete a summary report table.		
Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Does this standard include acceptance criteria? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
If no, include the results of testing in the 510(k).		
Does this standard include more than one option or selection of tests? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, report options selected in the summary report table.		
Were there any deviations or adaptations made in the use of the standard? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) <sup>5</sup> ?		
	<input type="checkbox"/>	<input type="checkbox"/>
Were deviations or adaptations made beyond what is specified in the FDA SIS? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, report these deviations or adaptations in the summary report table.		
Were there any exclusions from the standard? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, report these exclusions in the summary report table.		
Is there an FDA guidance <sup>6</sup> that is associated with this standard? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
If yes, was the guidance document followed in preparation of this 510k? .....		
	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Title of guidance: <u>Use of International Standard ISO-10993, 'Biological Evaluation of Medical Devices Part 1: Evaluation and Testing'</u> (Replaces #G87-1 #8294) (blue book memo), May 1, 1995		

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to his standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

## Traditional 510(k) Submission

### Straumann® Variobase™ Abutments

#### Declarations of Conformity and Summary Reports

(b) (4)		(b)(4) Trade Secret Process - Product Specs
A.	Applicable recognized consensus standard:	
B.	Requirements met?	Yes
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change
D.	Requirements that were not applicable to the device:	N/A
E.	Deviations from each applicable standard that were met including justification:	N/A
F.	Differences exist, if any, between the tested device to be marketed:	No
G.	Test laboratory:	(b)(4) Trade Secret Process - Product Specs

Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION  
 Traditional       Special       Abbreviated

STANDARD TITLE<sup>1</sup>  
**(b)(4) Trade Secret Process - Product Specs**

**Please answer the following questions** Yes      No

Is this standard recognized by FDA<sup>2</sup>? .....      

FDA Recognition number<sup>3</sup> ..... **#2-153**

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report<sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
 If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
 If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
 If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
 If yes, were deviations in accordance with the FDA supplemental information sheet (SIS)<sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
 If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
 If yes, report these exclusions in the summary report table.

Is there an FDA guidance<sup>6</sup> that is associated with this standard? .....         
 If yes, was the guidance document followed in preparation of this 510(k)? .....      

Title of guidance: Use of International Standard ISO-10993, 'Biological Evaluation of Medical Devices Part 1: Evaluation and Testing' (Replaces #G87-1 #8294) (blue book memo), May 1, 1995

<p><sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]</p> <p><sup>2</sup> Authority [21 U.S.C. 360d], <a href="http://www.fda.gov/cdrh/stdsprog.html">www.fda.gov/cdrh/stdsprog.html</a></p> <p><sup>3</sup> <a href="http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm">http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm</a></p> <p><sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or</p>	<p>certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.</p> <p><sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <a href="http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm">http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm</a></p> <p><sup>6</sup> The online search for CDRH Guidance Documents can be found at <a href="http://www.fda.gov/cdrh/guidance.html">www.fda.gov/cdrh/guidance.html</a></p>
---	--

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Declarations of Conformity and Summary Reports

(b)(4) Trade Secret Process - Product Specs

A.	Applicable recognized consensus standard:	(b)(4) Trade Secret Process - Product Specs
B.	Requirements met?	Yes
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change
D.	Requirements that were not applicable to the device:	N/A
E.	Deviations from each applicable standard that were met including justification:	N/A
F.	Differences exist, if any, between the tested device to be marketed:	No
G.	Test laboratory:	(b)(4) Trade Secret Process - Product Specs

Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional                       Special                       Abbreviated

STANDARD TITLE<sup>1</sup>

**(b)(4) Trade Secret Process - Product Specs**

**Please answer the following questions**

Yes      No

Is this standard recognized by FDA<sup>2</sup>? .....      

FDA Recognition number<sup>3</sup> ..... **#2-135**

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report<sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
 If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
 If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
 If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
 If yes, were deviations in accordance with the FDA supplemental information sheet (SIS)<sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
 If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
 If yes, report these exclusions in the summary report table.

Is there an FDA guidance<sup>6</sup> that is associated with this standard? .....         
 If yes, was the guidance document followed in preparation of this 510k? .....      

Title of guidance: Use of International Standard ISO-10993, 'Biological Evaluation of Medical Devices Part 1: Evaluation and Testing' (Replaces #G87-1 #8294) (blue book memo), May 1, 1995

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Declarations of Conformity and Summary Reports

(b)(4) Trade Secret Process - Product Specs		(b)(4) Trade Secret Process - Product Specs	
A.	Applicable recognized consensus standard:		
B.	Requirements met?	Yes	
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change	
D.	Requirements that were not applicable to the device:	N/A	
E.	Deviations from each applicable standard that were met including justification:	N/A	
F.	Differences exist, if any, between the tested device to be marketed:	No	
		(b)(4) Trade Secret Process - Product Specs	
G.	Test laboratory:		

Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional                       Special                       Abbreviated

STANDARD TITLE<sup>1</sup>

**(b)(4) Trade Secret Process - Product Specs**

**Please answer the following questions**

Yes      No

Is this standard recognized by FDA<sup>2</sup>? .....      

FDA Recognition number<sup>3</sup> ..... # \_\_\_\_\_

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report<sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
 If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
 If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
 If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
 If yes, were deviations in accordance with the FDA supplemental information sheet (SIS)<sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
 If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
 If yes, report these exclusions in the summary report table.

Is there an FDA guidance<sup>6</sup> that is associated with this standard? .....         
 If yes, was the guidance document followed in preparation of this 510k? .....      

Title of guidance: Use of International Standard ISO-10993, 'Biological Evaluation of Medical Devices Part 1: Evaluation and Testing' (Replaces #G87-1 #8294) (blue book memo), May 1, 1995

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)



# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Declarations of Conformity and Summary Reports

(b)(4) Trade Secret Process  
Product Specs

A.	Applicable recognized consensus standard:	(b)(4) Trade Secret Process - Product Specs
B.	Requirements met?	Yes
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change
D.	Requirements that were not applicable to the device:	N/A
E.	Deviations from each applicable standard that were met including justification:	N/A
F.	Differences exist, if any, between the tested device to be marketed:	No
G.	Test laboratory:	(b)(4) Trade Secret Process - Product Specs

Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION  
 Traditional       Special       Abbreviated

STANDARD TITLE<sup>1</sup>  
**(b)(4) Trade Secret Process - Product Specs**

Please answer the following questions	Yes	No
Is this standard recognized by FDA <sup>2</sup> ? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
FDA Recognition number <sup>3</sup> .....	<b>#8-63</b>	
Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Is a summary report <sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
If no, complete a summary report table.		
Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Does this standard include acceptance criteria? .....	<input checked="" type="checkbox"/>	<input type="checkbox"/>
If no, include the results of testing in the 510(k).		
Does this standard include more than one option or selection of tests? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, report options selected in the summary report table.		
Were there any deviations or adaptations made in the use of the standard? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) <sup>5</sup> ? .....		
Were deviations or adaptations made beyond what is specified in the FDA SIS? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, report these deviations or adaptations in the summary report table.		
Were there any exclusions from the standard? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, report these exclusions in the summary report table.		
Is there an FDA guidance <sup>6</sup> that is associated with this standard? .....	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, was the guidance document followed in preparation of this 510k? .....		
Title of guidance: _____		

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Declarations of Conformity and Summary Reports

#### 9.6 ISO 17664:2004

A.	Applicable recognized consensus standard:	(b)(4) Trade Secret Process - Product Specs
B.	Requirements met?	Yes
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change
D.	Requirements that were not applicable to the device:	N/A
E.	Deviations from each applicable standard that were met including justification:	N/A
F.	Differences exist, if any, between the tested device to be marketed:	No
G.	Test laboratory:	(b)(4) Trade Secret Process - Product Specs

Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional                       Special                       Abbreviated

STANDARD TITLE<sup>1</sup>

**(b)(4) Trade Secret Process - Product Specs**

**Please answer the following questions**

Yes      No

Is this standard recognized by FDA<sup>2</sup>? .....      

FDA Recognition number<sup>3</sup> ..... # \_\_\_\_\_

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report<sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
 If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
 If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
 If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
 If yes, were deviations in accordance with the FDA supplemental information sheet (SIS)<sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
 If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
 If yes, report these exclusions in the summary report table.

Is there an FDA guidance<sup>6</sup> that is associated with this standard? .....         
 If yes, was the guidance document followed in preparation of this 510k? .....         
 Title of guidance: \_\_\_\_\_

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Declarations of Conformity and Summary Reports

(b)(4) Trade Secret Process  
Product Specs

A.	Applicable recognized consensus standard:	(b)(4) Trade Secret Process - Product Specs
B.	Requirements met?	Yes
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change
D.	Requirements that were not applicable to the device:	N/A
E.	Deviations from each applicable standard that were met including justification:	N/A
F.	Differences exist, if any, between the tested device to be marketed:	No
G.	Test laboratory:	(b)(4) Trade Secret Process - Product Specs

Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional                       Special                       Abbreviated

STANDARD TITLE<sup>1</sup>

(b)(4) Trade Secret Process - Product Specs

**Please answer the following questions**

Yes      No

Is this standard recognized by FDA<sup>2</sup>? .....      

FDA Recognition number<sup>3</sup> ..... #14-261

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report<sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
 If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
 If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
 If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
 If yes, were deviations in accordance with the FDA supplemental information sheet (SIS)<sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
 If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
 If yes, report these exclusions in the summary report table.

Is there an FDA guidance<sup>6</sup> that is associated with this standard? .....         
 If yes, was the guidance document followed in preparation of this 510k? .....         
 Title of guidance: \_\_\_\_\_

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Declarations of Conformity and Summary Reports

(b)(4) Trade Secret Process - Product Specs

A.	Applicable recognized consensus standard:	(b)(4) Trade Secret Process - Product Specs
B.	Requirements met?	Yes
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change
D.	Requirements that were not applicable to the device:	N/A
E.	Deviations from each applicable standard that were met including justification:	N/A
F.	Differences exist, if any, between the tested device to be marketed:	No
G.	Test laboratory:	(b)(4) Trade Secret Process - Product Specs

(b)(4) Trade Secret Process - Product Specs

Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional                       Special                       Abbreviated

STANDARD TITLE<sup>1</sup>

**(b)(4) Trade Secret Process - Product Specs**

**Please answer the following questions**

Yes      No

Is this standard recognized by FDA<sup>2</sup>? .....      

FDA Recognition number<sup>3</sup> ..... # \_\_\_\_\_

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report<sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
 If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
 If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
 If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
 If yes, were deviations in accordance with the FDA supplemental information sheet (SIS)<sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
 If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
 If yes, report these exclusions in the summary report table.

Is there an FDA guidance<sup>6</sup> that is associated with this standard? .....         
 If yes, was the guidance document followed in preparation of this 510k? .....         
 Title of guidance: \_\_\_\_\_

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)



## Traditional 510(k) Submission

### Straumann® Variobase™ Abutments

#### Declarations of Conformity and Summary Reports

(b)(4) Trade Secret  
Process - Product Specs

A.	Applicable recognized consensus standard:	<b>ISO 14801:2007, Dentistry – Implants – Dynamic fatigue test for endosseous dental implants</b>
B.	Requirements met?	Yes
C.	Way(s) in which the standard may have been adapted (such as change in test method):	No change
D.	Requirements that were not applicable to the device:	N/A
E.	Deviations from each applicable standard that were met including justification:	N/A
F.	Differences exist, if any, between the tested device to be marketed:	No
G.	Test laboratory:	(b)(4) Trade Secret Process - Product Specs

Department of Health and Human Services  
 Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
 (To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION  
 Traditional       Special       Abbreviated

STANDARD TITLE<sup>1</sup>  
**(b)(4) Trade Secret Process - Product Specs**

**Please answer the following questions** Yes      No

Is this standard recognized by FDA<sup>2</sup>? .....      

FDA Recognition number<sup>3</sup> ..... **#4-195**

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report<sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
 If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
 If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
 If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
 If yes, were deviations in accordance with the FDA supplemental information sheet (SIS)<sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
 If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
 If yes, report these exclusions in the summary report table.

Is there an FDA guidance<sup>6</sup> that is associated with this standard? .....         
 If yes, was the guidance document followed in preparation of this 510(k)? .....      

Title of guidance: Guidance for Industry and FDA Staff 'Class II Special Controls Guidance Document: Root-form Endosseous Dental Implants and Endosseous Dental Implant Abutments', May 12, 2004

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Device Description

---

## 10 Device Description

The modification proposed in this premarket notification is to the Indications for Use for the previously cleared Straumann® CARES® Variobase™ Abutments for NNC, RN, WN, NC, and RC (K120822). The modified Indications for Use would allow Straumann to market the Straumann® Variobase™ Abutment as a stand-alone component. The dental laboratory would then manufacture the respective coping (the second component of the Variobase two-piece abutment) and/or prosthetic restoration via their preferred workflow of pressing, casting, or in-lab milling using a burnout coping or STL model for open CAD software.

### 10.1 Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Device Description

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Device Description

(b)(4) Trade Secret Process - Product Specs



**10.2 Basal Screws**

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Device Description

(b)(4) Trade Secret Process - Product Specs



**10.3 Accessories**

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Device Description

(b)(4) Trade Secret Process - Product Specs



**10.4 Procedure**

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Device Description

(b)(4) Trade Secret Process - Product Specs





**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Substantial Equivalence Discussion

---

## **11 Substantial Equivalence Discussion**

Within the meaning of the Medical Device Amendments Act of 1976, the proposed change to the Indications for Use for the Straumann® Variobase™ Abutments in this 510(k) premarket notification are substantially equivalent to the medical devices currently in commercial distribution listed below:

- K120822, Straumann CARES Variobase Abutment NNC, RN, WN, NC, RC
- K111935, Ti-Base Abutment (NT-Trading GmbH & Co. KG)

Straumann® Variobase™ Abutments are pre-manufactured (stock) abutments intended to be placed onto Straumann dental implants to provide support for customized prosthetic reconstructions such as crowns and bridges. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.

The NT-Trading Ti-Base Abutment is a pre-manufactured abutment supplied in two parts, the abutment and screw, for fixation onto dedicated endosseous dental implants and is intended for use as an aid in prosthetic restoration. As with the subject device, the coping/restoration is provided by the dental laboratory. The Ti-Base is compatible with commercially available dental CAD/CAM systems, such as 3Shape, Exocad, and Dental Wings. Such systems must be validated by the user.

The NT-Trading catalog and package insert are included in this submission in Appendices 3 and 4, respectively. The products that were cleared in premarket notification K111935 are outlined in Table 7 (the 510(k) Summary is included in Appendix 5). Specifically, the L-Series and N-Series abutments are compatible with implants of the Straumann Dental Implant System as shown in Table 7.

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Substantial Equivalence Discussion

---

(b)(4) Trade Secret Process - Product Specs



## Traditional 510(k) Submission

### Straumann® Variobase™ Abutments

#### Substantial Equivalence Discussion

The table below provides a comparison matrix of the proposed and predicate devices (K120822):

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
<b>K Number</b>	<b>Straumann® Variobase™ Abutments Subject Submission</b>	<b>Straumann® CARES® Variobase™ Abutments (K120822)</b>	
<b>Indications for Use</b>	<p>The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.</p>	<p>The Straumann® CARES® Variobase™ Abutment is a two-piece dental abutment consisting of the Straumann® Variobase™ Abutment and the Straumann® CARES® Variobase™ Coping which is intended to be placed onto Straumann dental implants to provide support for prosthetic reconstruction such as crown and bridges. Straumann® CARES® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.</p> <p>The Straumann® CARES® Variobase™ Coping polycon® ae in combination with the Straumann® CARES® Variobase™ Abutment is indicated for temporary (up to 180 days) dental restoration of a Straumann dental implant.</p>	<p>The indications for use are being modified to allow Straumann to market the Straumann® Variobase™ Abutment as a stand-alone component. The dental laboratory would then manufacture the respective coping and/or prosthetic restoration using a burnout coping or STL model for open CAD software.</p>
<b>Material</b>	Titanium-Aluminum-Niobium alloy (Ti-6Al-7Nb)	Titanium-Aluminum-Niobium alloy (Ti-6Al-7Nb)	Identical
<b>Abutment Diameter</b>	(b)(4) Trade Secret Process - Product Specs	(b)(4) Trade Secret	Identical
<b>Abutment Height</b>	(b)(4) Trade Secret Process - Product Specs	(b)(4) Trade Secret	Identical

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Substantial Equivalence Discussion

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
<b>K Number</b>	<b>Straumann® Variobase™ Abutments Subject Submission</b>	<b>Straumann® CARES® Variobase™ Abutments (K120822)</b>	
<b>Mode of Action</b>	Screw-retained or cement retained	Screw-retained or cement retained	Identical
<b>Reusable</b>	No	No	Identical

**Table 8 - Comparison Matrix: Proposed Device versus Predicate Devices (K120822)**

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Substantial Equivalence Discussion

The table below provides a comparison matrix of the proposed and predicate devices (K111935):

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
<b>K Number</b>	<b>Straumann® Variobase™ Abutments Subject Submission</b>	<b>Ti-Base Abutment (K111935)</b>	
<b>Indications for Use</b>	<p>The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.</p>	<p>Ti-Base for individual Zirconium Abutments: The devices covered by this submission are abutments which are placed into a dental implant to provide support for a prosthetic restoration.</p> <p>The Ti-Base for individual Zirconium Abutments are intended for use to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support single and multiple tooth prosthesis, in the mandible or maxilla. The prosthesis can be cement retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.</p> <p>The Ti-Base abutments are indicated for use with the following implant systems:</p> <ul style="list-style-type: none"> <li>• Nobel Biocare® Replace Select®</li> <li>• Nobel Biocare NobelActive™</li> <li>• Biomet 3i® Osseotite®</li> <li>• Biomet 3i® Osseotite® Certain®</li> <li>• Nobel Biocare Branemark®</li> <li>• Straumann® synOcta®</li> <li>• Straumann® Bone Level®</li> <li>• Zimmer® Tapered Screw-vent®</li> <li>• Astra Tech OsseoSpeed®</li> <li>• Dentsply-Friadent® Frialit®</li> </ul>	<p>Equivalent – Both the subject and predicate devices are designed to interface with the Straumann Bone Level or the Straumann Tissue Level implants.</p>

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Substantial Equivalence Discussion

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
<b>K Number</b>	<b>Straumann® Variobase™ Abutments Subject Submission</b>	<b>Ti-Base Abutment (K111935)</b>	
<b>Material</b>	Titanium-Aluminum-Niobium alloy (Ti-6Al-7Nb)	Titanium-Aluminum-Vanadium alloy (Ti-6Al-4V)	Equivalent
<b>Abutment Diameter</b>	(b)(4) Trade S t	(b)(4) Trade S t	Equivalent
<b>Abutment Height</b>	(b)(4) Trade S t	(b)(4) T d	Equivalent
<b>Mode of Action</b>	Screw-retained or cement retained	Screw-retained or cement retained	Identical
<b>Reusable</b>	No	No	Identical

**Table 9 - Comparison Matrix: Proposed Device versus Predicate Devices (K111935)**

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Proposed Labeling

---

## 12 Proposed Labeling

### 12.1 Package Label

There are no changes to the package label as a result of the proposed change in this premarket notification. To aid in the review of the submission, an example of the label is shown in Figure 1.





Figure 1 - Example label for Straumann® Variobase™ Abutment

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

Proposed Labeling

---

### 12.2 Proposed Package Insert/Instructions for Use for Straumann® Variobase™ Abutments

Instructions for use: Straumann® Variobase™ Abutments




Institut Straumann AG, Peter-Merian-Weg 12, CH-4002 Basel/Switzerland,

[www.straumann.com](http://www.straumann.com)

**English**

**CAUTION: Federal law restricts this device to sale by or on the order of a dental professional.**



#### 1. Product Description

##### Abutments

Abutments are placed into dental implants to provide support for prosthetic reconstructions such as crowns and bridges.

##### Basal Screws

Basal screws are used for the fixation of the abutment to the dental implant.

#### 2. Intended use

Prosthetic components directly connected to the endosseous dental implant are intended for use as an aid in prosthetic rehabilitations.





# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Proposed Labeling

The Straumann® Variobase™ Abutments have not been evaluated for safety and compatibility in the magnetic resonance environment. The Straumann® Variobase™ Abutments have not been tested for heating or migration in the magnetic resonance environment.

#### 6. Compatibility information

Straumann® implants and the prosthetic components are available in a variety of configurations to meet your clinical needs. The label on each product uses abbreviations to help you identify whether a particular abutment or coping is compatible with the implant that you are restoring. The implant as well as the prosthetic component contains an identifier for the connection type, summarized in the table below.

Implant connection type	Compatible prostheses
NC (Narrow CrossFit®)	parts labeled NC
RC (Regular CrossFit®)	parts labeled RC
NNC (Narrow Neck CrossFit®)	parts labeled NNC
RN (Regular Neck)	parts labeled RN
WN (Wide Neck)	parts labeled WN

#### 7. Cleaning and Disinfection

Straumann® Variobase™ Abutments and components are non-sterile when delivered. Before placing the restoration in the patient's mouth, the product must be cleaned, disinfected and sterilized. Straumann recommends the following procedure for cleaning, disinfection and sterilization of abutments prior to use.

- 1) Clean rinsing under flowing water while brushing outer and inner side with adequate brushes.
- 2) The pre-treated product is to be cleaned/disinfected in an automated washer disinfector. Select the appropriate program according to the manufacturer's instructions.

#### 8. Sterilization

The restoration may be sterilized unwrapped or can be placed in an accessory cassette and packaged twice in common sterilization wraps (paper/film bags). Steam sterilize according to the parameters below:

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Proposed Labeling

Material	Method	Conditions
Variobase™ Abutment, TAN Screw, TAN	Autoclave (moist heat) Displacement: gravity or fractionated vacuum	134 °C (273 °F) 5 minutes

**Please note:** User should ensure the use of the appropriate biological indicator for the sterilizer and parameters used.

**Please note:** User should consult the coping/restoration material manufacturer's recommendations regarding sterilization.

**Caution:** Use devices immediately after sterilization. Do not store sterilized devices.

### 9. Procedure

#### Use and handling of the Straumann® Variobase™ Abutments for the Dental Technician

##### Restoration design

Make a coping or crown following standard procedure according to the material manufacturer's instructions. When using pressing or casting techniques via wax-up, use the burn-out coping for Variobase™ Abutments which supports a clean and sharp-edged finish of the screw channel and a good fit to the Straumann® Variobase™ Abutment.

When using a digital workflow, use the Straumann Variobase Implant Kit with any software platform, to facilitate the precise design of the interface between the Straumann® Variobase™ Abutment and the coping. The kit consists of an STL file containing the required milling template for the inner coping geometry.

**Please note:** The following framework wall thickness guidelines must be followed:

For ceramic materials a framework wall thickness  $\geq 0.4$  mm

For polymer materials a framework wall thickness  $\geq 0.5$  mm

**Please note:** The materials used to fabricate the coping must have a flexural strength between 91.5 MPa and 1007 MPa.

##### Processing

Process the coping or crown following standard procedure according to the material manufacturer's instructions. Always finalize the crown or coping prior to bonding to the Straumann® Variobase™ Abutments.

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Proposed Labeling

---

#### Bonding

**Please note:** It is not necessary to sandblast the Straumann® Variobase™ Abutment.

- 1) Fix the abutment to the implant analog with a screw (hand-tight).
- 2) Seal the screw channel with wax.
- 3) Apply self-adhesive dental cement on the abutment. Only suitable self-adhesive cementation systems for the material used shall be used. Follow the instructions for use of both the dental material and cement/bonding material manufacturer. (Straumann® recommends Panavia™ F2.0 resin cement by Kuraray)
- 4) Bond the coping to the abutment.
- 5) Immediately remove excess cement from the abutment and polish the lower margin of the coping after the cement is set.
- 6) Optional: For cement retained-restorations: Make a crown following standard procedure according to the material manufacturer's instructions and finalize it.
- 7) Clean the restoration prior to sending it to the dentist.
- 8) Include this instruction for use when sending the restoration to the dentist.

#### **Use and handling of the Straumann® Variobase™ Abutments for the Dentist**

##### Remove the restoration from the master cast or the analog.

Clean, disinfect and sterilize the device as described in sections 7 and 8 of this Instructions for Use document.

##### Placing the restoration

- a) Remove the healing cap or temporary restoration.
- b) Clean and dry the interior of the implant and the abutment thoroughly.
- c) Place the sterilized restoration into the patient's mouth.
- d) Make sure that the retentive elements of the implant abutment connection are properly aligned.
- e) Use the screw delivered with the abutment to screw the abutment into the dental implant.

**Please note:** Always ensure that the surfaces of threads and screw heads are clean and that a new screw is used for the restoration.

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Proposed Labeling

f) Straumann® abutments are fixed to the implant using the Straumann® SCS screwdriver, ratchet and torque control device. Use the respective torque according to the table below:

Device type	Tightening torque	Special considerations
Abutments (permanent)	35 Ncm	n/a
Temporary abutments	15 – 35 Ncm	Tighten only to 35 Ncm if the implant is fully osseointegrated
Components on implant analogs	Hand-tight	n/a

For cement-retained restorations (optional):

- g) Close the screw channel with cotton and sealing compound (i.e., gutta-percha)
- h) Apply self-adhesive dental cement on the two-piece abutment. Only suitable self-adhesive cementation systems for the used materials shall be used. Follow the instructions for use of the cement/bonding material manufacturer (Straumann recommends Panavia™ F2.0 resin cement by Kuraray).
- i) Bond the crown to the two-piece abutment.
- j) Immediately remove excess cement from the two-piece abutment.

### Warning

Torques greater than 35 Ncm may result in failure of the abutment and/or implant. Torque values less than the recommended values may result in loosening of the abutment, which may lead to abutment and/or implant failure.

### Please note

Once the Straumann® abutment has been secured to the implant using the indicated torque, it should not be removed.

### 10. Further Information

For additional information about the use of Straumann® products, call Straumann's customer service department or visit [www.straumann.com](http://www.straumann.com).

For additional information, consult:

Basic information on the Straumann® Variobase™ Abutment

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

Proposed Labeling

---

### 11. Please note

Practitioners must have appropriate knowledge and instruction in the handling of the Straumann product described herein (“Straumann Product”) for using the Straumann Product safely and properly in accordance with these instructions for use.

The Straumann Product must be used in accordance with the instructions for use provided by the manufacturer. It is the practitioner’s responsibility to use the device in accordance with these instructions for use and to determine, if the device fits to the individual patient situation.

The Straumann Product is part of an overall concept and must be used only in conjunction with the corresponding original components and instruments distributed by Institut Straumann AG, its ultimate parent company and all affiliates or subsidiaries of such parent company (“Straumann”), except if stated otherwise in these instructions for use. If use of products made by third parties is not recommended by Straumann in these instructions for use, any such use will void any warranty or other obligation, express or implied, of Straumann.

### 12. Validity

Upon publication of these instructions for use, all previous versions are superseded.

© Institut Straumann AG, 2013. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

Panavia™ is a trademark of Kurary Co, LTD, JP.

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

Proposed Labeling

---

### Availability

Some items of the Straumann® Dental Implant System are not available in all countries.



Straumann Products with the CE mark fulfill the requirements of the Medical Devices Directive 93/42 EEC /



Refer to package insert



Manufacturer



Article number



Lot Number



Do not re-use



Non-sterile

Rx only

Federal law restricts this device to sale by or on the order of a dental professional.

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
**Sterilization and Shelf Life**

---

**13 Sterilization and Shelf Life**

**13.1 Sterilization**

(b)(4) Trade Secret Process - Product Specs





**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
**Sterilization and Shelf Life**

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
**Sterilization and Shelf Life**

---

**13.2 Shelf Life**

(b)(4) Trade Secret Process - Product Specs



**13.3 Packaging**

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Biocompatibility

---

**14 Biocompatibility**

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Biocompatibility

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Software

---

**15 Software**

Software requirements do not apply to Straumann® Variobase™ Abutments and basal screws because these are standard stock products that do not contain software. The elements of the Straumann predicate device (K120822) that are designed using the referenced CAD software are the zirconia copings. These copings are not included in the subject submission. Therefore, this section is not applicable as the proposed change does not contain and is not dependent on the use of software.

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

Electromagnetic Compatibility and Electrical Safety

---

### 16 Electromagnetic Compatibility and Electrical Safety

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Performance Testing – Bench

---

**17 Performance Testing – Bench**

(b)(4) Trade Secret Process - Testing



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Performance Testing – Bench

(b)(4) Trade Secret Process - Testing





**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Performance Testing – Bench

(b)(4) Trade Secret Process - Testing



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Performance Testing – Bench

(b)(4) Trade Secret Process - Testing



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Performance Testing – Bench

(b)(4) Trade Secret Process - Testing



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Performance Testing – Bench

(b)(4) Trade Secret Process - Testing



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Performance Testing – Bench

(b)(4) Trade Secret Process - Testing



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Performance Testing – Bench

(b)(4) Trade Secret Process - Testing



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Performance Testing – Bench

(b)(4) Trade Secret Process - Testing



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Performance Testing – Bench

(b)(4) Trade Secret Process - Testing





**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Performance Testing – Bench

(b)(4) Trade Secret Process - Testing



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Performance Testing – Animal

---

**18 Performance Testing – Animal**

This section is not applicable as no animal testing was performed in the development of the proposed devices.

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Performance Testing – Clinical

---

**19 Performance Testing – Clinical**

This section is not applicable as clinical study results are not being submitted in this premarket notification.

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 1

---

**Appendix 1 – Engineering Drawings**




























































**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 3


---

**Appendix 3 – NT-Trading Product Catalog**



# MORE THAN ABUTMENT'S

## CATALOGUE 2012/1





# CONTENTS

## MORE THAN ABUTMENT OPTIMIZED SAFETY AND QUALITY CONSTANT LOAD ACCORDING TO DIN EN ISO 14801

To analyze the behavior of the material and the design of the implants on so many of these force effects and to determine the load limits, they will be reproduced in simulations. The international standard ISO 14801 describes the test setup and the implementation of such a test of fatigue limit endosseous dental implants.


We have had constant load tests for endosseous dental implants carried out for n-trading abutments by the National Institute of Research and Development and Measurement Technique in accordance with standard DIN EN ISO 14801. The load corresponded to 30° angle to the implant axis. In addition, we have, in conformity with DIN EN ISO 14801 had a finite element analysis test carried out.

**NATIONAL INSTITUTE OF RESEARCH AND DEVELOPMENT FOR MECHATRONICS AND MEASUREMENT TECHNIQUE**

**TESTING REPORT**  
No. 113 from 24.06.2010

**APPROVED, GENERAL MANAGER**  
Univ.-Prof. PhD. eng. Gh. Ion GHEORGHE

F.P.OS - B.2.4 - 2



### IMPLANT SYSTEMS

<b>E-Serie</b> is compatible to Nobel Biocare Replace Select®	4
<b>F-Serie</b> is compatible to Nobel Biocare™ Nobel Active™	9
<b>H-Serie</b> is compatible to Biomet 3i Certain®	11
<b>I-Serie</b> is compatible to Biomet 3i Osseotite®	13
<b>K-Serie</b> is compatible to Nobel Biocare Brånemark®	14
<b>L-Serie</b> is compatible to Straumann Bone Level®	15
<b>N-Serie</b> is compatible to Straumann SynOcta®	19
<b>R-Serie</b> is compatible to Zimmer Tapered Screw-Vent®	24
<b>S-Serie</b> is compatible to Astra Tech OsseoSpeed®	28
<b>T-Serie</b> is compatible to Dentsply-Friadent Frialit/Xive®	31
<b>Prosthetic tools</b>	35
<b>Dentokeep</b>	36
<b>Scan equipment</b>	37
<b>nt-milling tools and grinding tools</b>	38
<b>Sinter components</b>	39
<b>Tightening torques, Conditions of warranty</b>	40
<b>Notes</b>	42
<b>NT-Preform</b>	43

nt-trading GmbH & Co. KG meets the requirements of  
ISO 13485:2007 and the EU Directive 93/42/EEC  
Health Canada Recognized Registrar CMCDCAS.



## E-SERIE

COMPATIBLE TO NOBEL BIOACARE REPLACE SELECT®




Platform  
3,5 mm NP

Platform  
4,3 mm RP

Platform  
5,0 mm WP

Platform  
6,0 mm

**Titanium Base**  
for individual milled  
Zirconium Abutment  
incl. screw  
Titan Grade 5



E 800

E 810

E 820

E 830

**Scan Body**  
for Titanium Base  
PEEK



E 00 W

E 10 W

E 20 W

E 30 W

**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK



E 9.3D3.500

E 9.3D4.300

E 9.3D5.000

E 9.3D6.000

**Lab Analog**  
Stainless Steel



E 50

E 51

E 52

E 53

**Abutment Screw**  
Titan Grade 5  
Recommended tightening  
torques 35 Ncm



E 60

E 61

E 61

E 61

## E-SERIE

COMPATIBLE TO NOBEL BIOACARE REPLACE SELECT®



Platform  
4,3 mm RP

Platform  
5,0 mm WP


**2-CONnect-Base Set**  
incl. 2-CONnect Cap and  
Cap-Screw  
Titan Grade 5



E 810 S

E 820 S

**2-CONnect-Base**  
Titan Grade 5  
Recommended tightening  
torques 35 Ncm



E 810 M

E 820 M

**2-CONnect Cap**  
Titan Grade 5



E 810 F

E 820 F

**Cap-Screw**  
Titan Grade 5  
Recommended tightening  
torques 30 Ncm



N 60

N 60

**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK



E 9.3D4.300

E 9.3D5.000

**Lab Analog**  
Stainless Steel



E 51

E 52

## E-SERIE

COMPATIBLE TO NOBEL BIOCARE REPLACE SELECT®

- 

Platform  
3,5 mm NP
- 

Platform  
4,3 mm RP
- 

Platform  
5,0 mm WP

**Straight Abutment**  
incl. Screw  
GH 1,0 mm



**Straight Abutment**  
incl. Screw  
GH 2,5 mm



## E-SERIE

COMPATIBLE TO NOBEL BIOCARE REPLACE SELECT®

- 

Platform  
3,5 mm NP
- 

Platform  
4,3 mm RP
- 

Platform  
5,0 mm WP

**Angled Abutment 16°**  
angled over surface  
GH 1,0 mm  
incl. Screw



**Angled Abutment 16°**  
angled over edge  
GH 1,0 mm  
incl. Screw



**Angled Abutment 16°**  
angled over surface  
GH 2,5 mm  
incl. Screw



**Angled Abutment 16°**  
angled over edge  
GH 2,5 mm  
incl. Screw





# E-SERIE

COMPATIBLE TO NOBEL BIOACARE REPLACE SELECT®



Platform  
3,5 mm NP

Platform  
4,3 mm RP

Platform  
5,0 mm WP



**HSL Abutment**  
notation indexed  
incl. Screw



**Implant Pic-Up**  
incl. screw

# F-SERIE

COMPATIBLE TO NOBEL BIOACARE™ NOBEL ACTIVE™



Platform  
3,5 mm NP

Platform  
4,3 mm / 5,0 mm RP



**Titanium Base**  
for individual milled  
Zirconium Abutment  
incl. screw  
Titan Grade 5



**Scan Body**  
for Titanium Base  
PEEK



**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK



**Lab Analog**  
Stainless Steel




**Abutment Screw**  
Titan Grade 5  
Recommended tightening  
torques 35 Ncm

# F-SERIE

KOMPATIBEL ZU NOBEL BIO CARE™ NOBEL ACTIVE™

Traditional 510(k)  
Straumann® Variobase™ Abutments



Platform  
3,5 mm NP

Platform  
4,3 mm / 5,0 mm RP

**2-CONnect-Base Set**  
incl. 2-CONnect Cap and Cap-Screw  
Titan Grade 5

	F 800 S
	F 810 S

**2-CONnect-Base**  
Titan Grade 5  
Recommended tightening torques 35 Ncm

	F 800 M
	F 810 M

**2-CONnect Cap**  
Titan Grade 5

	F 800 F
	F 810 F

**Cap-Screw**  
Titan Grade 5  
Recommended tightening torques 30 Ncm

	N 60
--	------

**Scan Body 3D Guide**  
for Titanium Base + 2-CONnect PEEK

	F 9.3D3.500
	F 9.3D4.350

**Lab Analog**  
Stainless Steel

	F 50
	F 51

# H-SERIE

KOMPATIBEL ZU BIOMET 3I CERTAIN®



Platform  
3,4 mm

Platform  
4,1 mm

Platform  
5,0 mm

**Titanium Base**  
for individual milled Zirconium Abutment incl. screw  
Titan Grade 5

	H 800
	H 810
	H 820

**Scan Body**  
for Titanium Base PEEK

	H 00 W
	H 10 W

**Scan Body 3D Guide**  
for Titanium Base + 2-CONnect PEEK

	H 9.3D3.400
	H 9.3D4.150
	H 9.3D4.150

**Lab Analog**  
Stainless Steel

	H 50
	H 51
	H 52

**Abutment Screw**  
Titan Grade 5  
Recommended tightening torques 20 Ncm

	H 60
	H 60
	H 60



# H-SERIE

COMPATIBLE TO BIOMET 3I CERTAIN®



Platform  
3,4 mm

Platform  
4,1 mm

<b>2-CONnect-Base Set</b> incl. 2-CONnect Cap and Cap-Screw Titan Grade 5		
	H 800 S	H 810 S

<b>2-CONnect-Base</b> Titan Grade 5 Recommended tightening torques 20 Ncm		
	H 800 M	H 810 M

<b>2-CONnect Cap</b> Titan Grade 5		
	H 800 F	H 810 F

<b>Cap-Screw</b> Titan Grade 5 Recommended tightening torques 15 Ncm		
	KS 60	KS 60

<b>Scan Body 3D Guide</b> for Titanium Base + 2-CONnect PEEK		
	H 9.3D3.400	H 9.3D4.150

<b>Lab Analog</b> Stainless Steel		
	H 50	H 51

# I-SERIE

COMPATIBLE TO BIOMET 3I OSSEOTITE®



Platform  
3,4 mm

Platform  
4,1 mm

Platform  
5,0 mm

<b>Titanium Base</b> for individual milled Zirconium Abutment incl. screw Titan Grade 5			
	I 800	I 810	I 820

<b>Scan Body</b> for Titanium Base PEEK			
	I 100 W	I 110 W	I 120 W

<b>Scan Body 3D Guide</b> for Titanium Base + 2-CONnect PEEK			
	I 9.3D3.400	I 9.3D4.150	I 9.3D4.150

<b>Lab Analog</b> Stainless Steel			
	I 50	I 51	I 52

<b>Straight Abutment</b> incl. Screw GH 2,5 mm			
	I 110	I 110	I 110

<b>Abutment Screw</b> Titan Grade 5 Recommended tightening torques 35 Ncm			
	I 61	I 61	I 61

## K-SERIE

COMPATIBLE TO NOBEL BIOCARE BRÄNEMARK®



Platform  
3,5 mm

Platform  
4,1 mm

Platform  
5,1 mm

**Titanium Base**  
for individual milled  
Zirconium Abutment  
incl. screw  
Titan Grade 5



K 800

K 810

K 820

**Scan Body**  
for Titanium Base  
PEEK



K 00 W

K 10 W

K 20 W

**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK



K 9.3D3.500

K 9.3D4.100

K 9.3D5.100

**Lab Analog**  
Stainless Steel



K 50

K 51

K 52

**Straight Abutment**  
incl. Screw  
GH 2,5 mm



K 110

**Abutment Screw**  
Titan Grade 5  
Recommended tightening  
torques 35 Ncm



K 60

K 61

K 62

## L-SERIE

COMPATIBLE TO STRAUMANN BONE LEVEL®



Platform  
3,3 mm NC

Platform  
4,1 mm / 4,8 mm RC

**Titanium Base**  
for individual milled  
Zirconium Abutment  
incl. screw  
Titan Grade 5



L 800

L 810

**Scan Body**  
for Titanium Base  
PEEK



L 00 W

L 10 W

**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK



L 9.3D3.300

L 9.3D4.148

**Lab Analog**  
Stainless Steel



L 50

L 51

**Abutment Screw**  
Titan Grade 5  
Recommended tightening  
torques 35 Ncm




L 60

L 60



# L-SERIE

COMPATIBLE TO STRAUMANN BONE LEVEL®



Platform  
3,3 mm NC

Platform  
4,1 mm / 4,8 mm RC

<b>2-CONnect-Base Set</b> incl. 2-CONnect Cap and Cap-Screw Titan Grade 5		
<b>2-CONnect-Base</b> Titan Grade 5 Recommended tightening torques 35 Ncm		
<b>2-CONnect Cap</b> Titan Grade 5		
<b>Cap-Screw</b> Titan Grade 5 Recommended tightening torques 30 Ncm		
<b>Scan Body 3D Guide</b> for Titanium Base + 2-CONnect PEEK		
<b>Lab Analog</b> Stainless Steel		

# L-SERIE

COMPATIBLE TO STRAUMANN BONE LEVEL®



Platform  
3,3 mm NC

Platform  
4,1 mm / 4,8 mm RC


<b>Straight Abutment</b> incl. Screw GH 2,5 mm		
<b>Straight Abutment</b> incl. Screw GH 3,0 mm		
<b>Angled Abutment 18°</b> angled over surface GH 1,5 mm incl. Screw		
<b>Angled Abutment 18°</b> angled over edge GH 1,5 mm incl. Screw		



# L-SERIE

COMPATIBLE TO STRAUMANN BONE LEVEL®

Traditional 510(k)  
Straumann® Variobase™ Abutments



Platform 3,3 mm NC  
Platform 4,1 mm / 4,8 mm PC

<b>Angled Abutment 18°</b> angled over surface GH 3,0 mm incl. Screw		L 200-1-3
		L 210-1-3

<b>Angled Abutment 18°</b> angled over edge GH 3,0 mm incl. Screw		L 200-2-3
		L 210-2-3

<b>HSL Abutment</b> rotation indexed incl. Screw		L 11.CA3.300
		L 11.CA4.148

<b>Implant Pic-Up</b> incl. screw		L TRAND013.3
		L TRAR0024.1

Notice: Products indicated with ® are registered brand names of the respective manufacturers.

# N-SERIE

COMPATIBLE TO STRAUMANN SYNOCTA®



Platform 3,5 mm NN  
Platform 4,8 mm RN  
Platform 6,5 mm WN

<b>Titanium Base</b> for individual milled Zirconium Abutment incl. screw Titan Grade 5		N 800
		N 810
		N 820

<b>Scan Body</b> for Titanium Base PEEK		N 00 W
		N 10 W
		N 20 W

<b>Scan Body 3D Guide</b> for Titanium Base + 2-CONnect PEEK		N 9.3D3.500
		N 9.3D4.800
		N 9.3D6.500

<b>Lab Analog</b> Stainless Steel		N 50
		N 51
		N 52

<b>Abutment Screw</b> Titan Grade 5 Recommended tightening torques 35 Ncm		N 60
		N 62
		N 62

Notice: Products indicated with ® are registered brand names of the respective manufacturers.

# N-SERIE

COMPATIBLE TO STRAUMANN SYNOCTA®



<b>2-CONnect-Base Set</b> incl. 2-CONnect Cap and Cap-Screw Titan Grade 5		
<b>2-CONnect-Base</b> Titan Grade 5 Recommended tightening torques 35 Ncm		
<b>2-CONnect Cap</b> Titan Grade 5		
<b>Cap-Screw</b> Titan Grade 5 Recommended tightening torques 30 Ncm		
<b>Scan Body 3D Guide</b> for Titanium Base + 2-CONnect PEEK		
<b>Lab Analog</b> Stainless Steel		

Notice: Products indicated with ® are registered brand names of the respective manufacturers.

# N-SERIE

COMPATIBLE TO STRAUMANN SYNOCTA®



<b>Straight Abutment</b> incl. Screw			
<b>Angled Abutment 16°</b> angled over surface incl. Screw			
<b>Angled Abutment 16°</b> angled over edge incl. Screw			
<b>Abutment Screw</b> Titan Grade 5 Recom. tight. torques 35 Ncm			
<b>RN-Massiv Abutment</b> only for Dentist			
<b>WN-Massiv Abutment</b> only for Dentist			

Notice: Products indicated with ® are registered brand names of the respective manufacturers.



## N-SERIE

COMPATIBLE TO STRAUMANN SYNOCTA®



Platform  
3,5 mm NN

Platform  
4,8 mm RN

Platform  
6,5 mm WN



Platform  
3,5 mm NN



Platform  
4,8 mm RN



Platform  
6,5 mm WN




## N-SERIE

COMPATIBLE TO STRAUMANN SYNOCTA®



**Straight Abutment**  
incl. Screw



**Angled Abutment 16°**  
incl. Screw


**Abutment Screw**



**HSL Abutment**  
notation indexed  
incl. Screw



**HSL Abutment**  
rotating  
incl. Screw




**Abutment Screw**  
Titan Grade 5  
Recommended tightening  
torques 35 Nm



**Angled Abutment 21°**  
angled over surface  
incl. Screw



**Angled Abutment 21°**  
angled over edge  
incl. Screw



**WN-Massiv Abutment**  
Height 4,0 mm



**Implant Pic-Up**  
incl. screw

# R-SERIE

COMPATIBLE TO ZIMMER TAPERED SCREW-VENT®



Platform  
3,5 mm

Platform  
4,5 mm

Platform  
5,7 mm

**Titanium Base**  
for individual milled  
Zirconium Abutment  
incl. screw  
Titan Grade 5



R 800

R 810

R 820

**Scan Body**  
for Titanium Base +  
PEEK



R 00 W

R 10 W

R 20 W

**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK



R 9.3D3.500

R 9.3D4.500

R 9.3D5.700

**Lab Analog**  
Stainless Steel



R 50

R 51

R 52

**Abutment Screw**  
Hex 0,50" (1,26 mm)  
Recommended tightening  
torques 30 Ncm



R 60

R 60

# R-SERIE

COMPATIBLE TO ZIMMER TAPERED SCREW-VENT®



Platform  
3,5 mm

Platform  
4,5 mm

**2-CONnect-Base Set**  
incl. 2-CONnect Cap and  
Cap-Screw  
Titan Grade 5



R 800 S

R 810 S

**2-CONnect-Base**  
Titan Grade 5  
Recommended tightening  
torques 30 Ncm



R 800 M

R 810 M

**2-CONnect Cap**  
Titan Grade 5



R 800 F

R 810 F

**Cap-Screw**  
Titan Grade 5  
Recommended tightening  
torques 25 Ncm



KS 60

KS 60

**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK



R 9.3D3.500

R 9.3D4.500

**Lab Analog**  
Stainless Steel



R 50

R 51

Notice: Products indicated with ® are registered brand names of the respective manufacturers.

Notice: Products indicated with ® are registered brand names of the respective manufacturers.





## R-SERIE

COMPATIBLE TO ZIMMER TAPERED SCREW-VENT®



Platform  
5,7 mm



Platform  
4,5 mm



Platform  
3,5 mm

### Straight Abutment

incl. Screw  
GH 2,5 mm



R 100



R 110



R 120

### Angled Abutment 16°

angled over surface  
incl. Screw  
GH 2,5 mm



R 200-1



R 210-1



R 220-1

### Angled Abutment 16°

angled over surface  
incl. Screw  
GH 2,5 mm



R 200-2



R 210-2



R 220-2

### Massiv Abutment

Titan Grade 5



R 400



R 410



R 420

### HSL Abutment

rotation indexed  
incl. Screw



R 300



R 310



R 320

### Abutment Screw

Hex 0,50" (1,26 mm)  
Recommended tightening  
torques 30 Ncm



R 60



R 60



R 60

Notice: Products indicated with ® are registered brand names of the respective manufacturers.

## R-SERIE

COMPATIBLE TO ZIMMER TAPERED SCREW-VENT®



Platform  
4,5 mm



Platform  
3,5 mm

### Implant Pic-Up

incl. screw



R TR-00013.5



R TR-00024.5





# S-SERIE

COMPATIBLE TO ASTRA TECH OSSEOSPEED®



Platform  
3,5 mm / 4,0 mm



Platform  
4,5 mm / 5,0 mm

**Titanium Base**  
for individual milled  
Zirconium Abutment  
incl. screw  
Titan Grade 5



**Scan Body**  
for Titanium Base  
PEEK



**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK



**Lab Analog**  
Stainless Steel



**Straight Abutment**  
incl. Screw  
GH 1,5 mm



**Abutment Screw**  
Titan Grade 5  
Recommended tightening  
torques 25 Ncm



# S-SERIE

COMPATIBLE TO ASTRA TECH OSSEOSPEED®



Platform  
3,5 mm / 4,0 mm



Platform  
4,5 mm / 5,0 mm

**2-CONnect-Base Set**  
incl. 2-CONnect Cap and  
Cap-Screw  
Titan Grade 5



**2-CONnect-Base**  
Titan Grade 5  
Recommended tightening  
torques 25 Ncm



**2-CONnect Cap**  
Titan Grade 5



**Cap-Screw**  
Titan Grade 5  
Recommended tightening  
torques 20 Ncm



**Scan Body 3D Guide**  
for Titanium Base +  
2-CONnect  
PEEK



**Lab Analog**  
Stainless Steel





# S

## S-SERIE

COMPATIBLE TO ASTRA TECH OSSEOSPEED®



Platform 3,5 mm / 4,0 mm    Platform 3,5 mm / 4,0 mm    Platform 4,5 mm / 5,0 mm

<b>Straight Abutment</b> incl. Screw GH 1,5 mm	S 110	S 120	
<b>Angled Abutment 16°</b> incl. Screw GH 1,5 mm	S 200	S 210	S 220
<b>HSL Abutment</b> notation indexed incl. Screw	S 11.CA3.540	S 11.CA4.550	
<b>Implant Pic-Up</b> incl. screw	S TR-00013.5	S TR-00024.5	

## T-SERIE

COMPATIBLE TO DENTSPLY-FRIADENT FRIALIT/XIVE®



Platform 3,4 mm    Platform 3,8 mm    Platform 4,5 mm    Platform 5,5 mm

<b>Titanium Base</b> for individual milled Zirconium Abutment incl. screw Titan Grade 5	T 800	T 805	T 810	T 820
<b>Scan Body</b> for Titanium Base PEEK	T 00W	T 05W	T 10W	T 10W
<b>Scan Body 3D Guide</b> for Titanium Base + 2-CONnect PEEK	T 9.3D3.400	T 9.3D3.800	T 9.3D4.555	T 9.3D4.555
<b>Lab Analog</b> Stainless Steel	T 50	T 55	T 51	T 52
<b>Abutment Screw</b> Titan Grade 5 Recommended tightening torques 25 Ncm	T 60	T 60	T 60	T 60

# T-SERIE

COMPATIBLE TO DENTSPLY-FRIADENT FRIALIT/XIVE®



<b>2-CONnect-Base Set</b> incl. 2-CONnect Cap and Cap-Screw Titan Grade 5			
<b>2-CONnect-Base</b> Titan Grade 5 Recommended tightening torques 25 Ncm			
<b>2-CONnect Cap</b> Titan Grade 5			
<b>Cap-Screw</b> Titan Grade 5 Recommended tightening torques 15 Ncm			
<b>Scan Body 3D Guide</b> for Titanium Base + 2-CONnect PEEK			
<b>Lab Analog</b> Stainless Steel			

Notice: Products indicated with ® are registered brand names of the respective manufacturers.

# T-SERIE

COMPATIBLE TO DENTSPLY-FRIADENT FRIALIT/XIVE®



<b>Angled Abutment 16°</b> angled over edge GH 1,0 mm		
<b>Angled Abutment 16°</b> angled over surface GH 1,0 mm		
<b>Angled Abutment 16°</b> angled over surface GH 2,5 mm		
<b>Straight Abutment</b> incl. Screw GH 1,0 mm		
<b>Straight Abutment</b> incl. Screw GH 2,5 mm		
<b>Abutment Screw</b> Hex 0.50" (1,26 mm) Recommended tightening torques 25 Ncm		

Notice: Products indicated with ® are registered brand names of the respective manufacturers.



# T-SERIE

COMPATIBLE TO DENTSPLY-FRIADENT FRIALIT/XIVE®



**Implant Pic-Up**  
incl. screw

T TR-00013.4

T TR-00023.8

T TR-00034.5



**HSL Abutment**  
notation indexed  
incl. Screw

T 11.CA3.400

T 11.CA3.800

T 11.CA4.500

# PROSTHETIC TOOLS



**Torque Ratchet**  
continuously adjustable  
to max. 40 Ncm  
ISO-Shaft Connection

W 11.000.000



**Laboratory**  
Screwdriver  
for changeable inserts  
(without ISO-Shaft)  
ISO-Shaft Connection

W 11.100.000



**Insert**  
for Laboratory Screwdriver  
and Torque Ratchet  
ISO-Shaft for S-, I- and R-Series  
Hex 0,50" (1,26 mm)

W 11.IRS.G10



**Insert**  
for Laboratory Screwdriver  
and Torque Ratchet  
ISO-Shaft for T- and H-Series  
Hex 1,20 mm

W 11.TH.O.G20



**Insert**  
for Laboratory Screwdriver  
and Torque Ratchet  
ISO-Shaft for Ankylos®  
Hex 1,00 mm

W 11.Y00.G30



**Insert**  
for Laboratory Screwdriver  
and Torque Ratchet  
ISO-Shaft for L- and N-Series  
Torx T6

W 11.LND.G40



**Insert**  
for Laboratory Screwdriver  
and Torque Ratchet  
ISO-Shaft for E-, F- and  
K-Series UG

W 11.EFK.G50



**Insert**  
for Laboratory Screwdriver  
and Torque Ratchet ISO-Shaft  
for 2-CONNECT-Base-Primary

W 11.005.G60

# DENTOKEEP PEEK DISC

## Dentokeep

nt-trading Dentokeep Disc PEEK are blanks for technical milling manufacture of permanent dentures and prosthetic brace in the CAM process. The derived prosthetic designs are available on the remaining teeth, - can be attached to stumps and / or implant abutment and provide a basic functional and aesthetic care.

1. PEEK Disc 98.5 mm, pearl white, 14 mm and 18 mm thick
2. For CAD / CAM Technology for CAD/CAM Technology
3. Clip in prosthetics CAD/CAM
4. Class II product for permanent prosthetic restorations



## Dentokeep

PEEK Disc  
Size 14 mm

12.000.P14



## Dentokeep

PEEK Disc  
Size 18 mm

12.000.P18



# SCAN EQUIPMENT

## nt-OptiScan™ Spray

nt-OptiScan™ Spray a new product development of nt-trading GmbH & Co. KG.

## nt-OptiScan™ Spray

### Advantages:

- Optimize composition
- Perfect dosage with proven applicator
- Optimize preparation of your working model for optical scanning with modern dental scanner
- Economical application
- Water soluble
- Handy size: Volume 75 ml brutto
- Superior scanning results
- Special attractive price





# NT-MILLING AND GRINDING TOOLS

nt-trading tools are specially developed for dental CAD/CAM systems. They are up to the high mark to process dental materials and are constant liable to quality checks.  
Tools coated with tax are especially qualified for process of dental high-capacity plastics, as well as for process of milling wax.  
Tools coated with diamond are especially qualified for process of abrasive materials like zirconium dioxide.

<b>nt-Diaburr</b>			
			
0,6 mm	1,0 mm	2,5 mm	
<b>11100</b>	<b>11101</b>	<b>11102</b>	

The milling tool coated with diamonds nt-Diaburr offers you a lifetime of at least 350 units while process of ZrO2 with the milling machine.

<b>nt-Uburr</b>			
			
1,0 mm	2,5 mm		
<b>11201</b>	<b>11202</b>		

Universal milling tools are qualified for process of plastics, waxes and ZrO2 with the milling machine. Depending on choice of material the lifetime is at least 175 units.

<b>nt-PEEKburr</b>			
			
1,0 mm	2,5 mm		
<b>11301</b>	<b>11302</b>		

Special milling tools are qualified for process of PEEK materials like Dentokeep, plastics and waxes.

# SINTER COMPONENTS

## nt-pearls

Sinterization of Zr-oxide substructures.

## Sinter pearls for dental CAD/CAM technology

### Advantages:

Sinterization of Zr-oxide substructures. Full ceramics are the answer to patients demands for highly esthetic, metal-free and durable restorations Zr-oxide became a common substructure for dental restorations already quite a while ago. This material features outstanding biocompatibility and mechanical properties Zr-oxide is mechanically machined in its greenstage, taking into consideration the specific shrinkage of this material. The sinterization is optimized for the specific material and its oversized milling relative to the material properties, thus assuring the best possible fit.

Another important factor contributing to a perfect fit is the right choice of support for this to be sintered substructure. Every manufacturer has his own recommendations concerning these supports. Some suggest to sinter on special „pearls“, unfortunately particle size and quality are not always right to by example avoid getting pinched between the interproximals and thus distorting the bridge.

Other solutions favour the use of „drops“, resting on very expensive and fragile support discs. The necessary preparations for this type of support are always time consuming, sometimes even cumbersome and allways generating additional costs.

A viable alternative are special hi-density and quality sinterization pearls from nt-trading, optimized for crowns and bridges and therefore avoiding the common „pinching“ and distortion problems especially known to bridges.





# INSTRUCTION FOR USE

## Indication:

For manufacturing of individual abutments on dental implants. The individual abutments can be combined with copings, crowns or suprastructures made of dental ceramics.

## Contraindication:

The Ti-Bases of each Series can only be combined with the matching implant, e.g. the E-Series shall be combined exclusively with Replace Select® Implants. They cannot be combined with implants of a different implant type or manufacturer. The diameter of the Ti-base must correspond in size to the used implant in order to prevent a peri-implant tissue irritation. The Ti-Bases are indicated for single use only. If they are used multiple times, they might damage the implants.

For fixation of the Ti-Bases on the implant, the correct torque force, recommended by the implant manufacturer, has to be considered carefully to avoid the damage of the implant-bone connection.

Ncm	Abutment
20	H-Series
25	S-Series T-Series F-Series
30	R-Series
35	I-Series K-Series N-Series E-Series L-Series

Mechanical treatment of the connection part of the Ti-Base will damage the correct fitting of the Ti-Base on the implant.

## Handling method:

**Ceramic abutments:**  
Milling with CAD/CAM-machines of zirconium oxide- or aluminum oxide- ceramics according to the anatomic form of a crown or coping. The ceramic copings or crowns shall be milled or polished with diamond instruments and with minimal pressure and water-cooling. The minimal thickness shall be 0.5 mm Sharpe edges must be avoided.

## Veneering:

Copings shall be veneered with appropriate ceramics before cementing onto the Ti-Base. The instructions for use of the ceramic manufacturers have to be considered. Treatment of the Ti-Base and the ceramic abutment before cementing: Sandblasting of the contact surfaces with Al<sub>2</sub>O<sub>3</sub>, 50 µm, 2 bar and intensive cleaning of dust and grease. It is recommended to protect the connection part of the Ti-Base with an implant analog during handling.

## Cementing:

It is recommended to cement the ceramic abutment onto the Ti-Base with Panavia® F2.0 (Kuraray) with RelayXUnicem® (3M Espe) or other equivalent cements. The instructions for use of the cements shall be followed carefully. The Ti-Base shall be fixed onto an implant analog with the abutments screw. The head of the screw has to be covered with wax or resin. The mixed cement is applied onto the contact part of the Ti-Base. The abutment is pressed onto the Ti-Base. The final position is evaluated by slight rotation. The gap between abutment and the Ti-Base must be as small as possible. Remaining cement shall be removed immediately.

## Polishing:

After hardening the remaining cement shall be removed with rotating silicon instruments. The cement inside the screw channel has to be removed carefully.

## Scan Body:

### Indication:

#### Scan Body:

For the CAD/CAM scanning of the model, the Scan Body is used to indicate the position of the implant. The size of the Scan Body shall be corresponding to the original Implant system, implant diameter and Ti-Base Series.

The chamfer of the Scan Body prevents the rotation of the ceramic abutment. The Scan Body is fixed on the implant analogue with the abutment screw. After correct positioning, there is no gap visible between implant and Scan Body. Rotation of the Scan Body is impossible

## Tightening torques

Ncm	2-CONnect-M-Abutment	
20	H-Series	
25	S-Series	T-Series
30	R-Series	
35	E-Series	F-Series L-Series N-Series

Ncm	2-CONnect Cap-Screw N 60	
20	S-Series	T 805, T 810
30	E-Series	F-Series L-Series N-Series

Ncm	2-CONnect Cap-Screw KS 60	
15	H 800, H 810	T 800 R 800 R 810

## Conditions of warranty

Within our general terms of sale and warranty we ensure the perfect quality of our products. Due to our high production standards can offer you a 10-years' warranty on our prosthetic components.

We offer you a warranty on our nitriding components according to the conditions of warranty.

The warranty includes all material and manufacturing defects which may occur within the 10-years' time of warranty. We only give warranty for our contracting/purchasing partners (dentists, dental hospitals, laboratory). Any other persons besides those mentioned cannot lay claims to the warranty. It is not possible to assign the warranty claims.





**The new address  
from February 2013:**

nt-trading GmbH & Co. KG  
Nördliche Uferstraße 8  
76189 Karlsruhe  
Germany

PCF/D-2012-12

nt-trading GmbH & Co. KG | Esso-Straße 16 | 76187 Karlsruhe | Germany  
Phone: +49 (0)721 / 91 54 71 - 60 | [info@nt-trading.com](mailto:info@nt-trading.com) | [www.nt-trading.com](http://www.nt-trading.com)

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 4

---

**Appendix 4 – NT-Trading Ti-Base Package Insert**

## Instruction for use

### Ti- Base-, for Individual Abutments, Scan Base and Scan Body E-, I-, K-, N-, R-, S-, T-, H-, L-, and F- Series

CE 0123

#### INDICATION:

The Ti-Base, is intended for use to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support single and multiple tooth prosthesis, in the mandible or maxilla. The prosthesis can be cement retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.

The Ti-Base are for manufacturing of individual abutments on dental implants. The individual abutments can be combined with restorations such as, crowns or superstructures made of dental ceramics.

The Ti-Bases of the E-Series are indicated for Replace select® implants, manufactured by Nobel Biocare®.

The Ti-Base of the I-Series are indicated for Osseotite® implants manufactured by Biomet 3i®.

The Ti-Bases of the K-Series are indicated for Brånemark® implants, manufactured by Nobel Biocare®.

The Ti-Bases of the N-Series are indicated for synOcta® implants, manufactured by Straumann®.

The Ti-Bases of the R-Series are indicated for Tapered Screw-vent® implants, manufactured by Zimmer®.

The Ti-Bases of the S-Series are indicated for OsseoSpeed® implants, manufactured by Astra Tech®.

The Ti-Bases of the T-Series are indicated for Frialit® implants, manufactured by Dentsply-Friadent®.

The Ti-Bases of the H-Series are indicated for Biomet 3i Certain® implants, manufactured by Biomet 3i®.

The Ti-Bases of the L-Series are indicated for Bone Level implants, manufactured by Straumann®.

The Ti-Bases of the F-Series are indicated for NobelActive™ implants, manufactured by Nobel Biocare®.

#### ARTICLE NO. OF TI-BASES AND ABUTMENT SCREWS

The Abutments Ti-Base E-Series is compatible with Nobel Biocare Replace Select® Implant.

Nobel Replace Select®	3.5 mm	4.3 mm	5.0 mm	6.0 mm
Ti-Base E-Series	E 800	E 810	E 820	E 830
Abutment screw	E 60	E 61	E 61	E 61

The Abutment Ti-Base I-Series is compatible with Biomet 3i Osseotite®.

Biomet 3i Osseotite®	3.4 mm	4.1 mm	5.0 mm
Ti-Base I-Series	I 800	I 810	I 820
Abutment screw	I 61	I 61	I 61

The Abutment Ti-Base K-Series is compatible with Nobel Biocare Brånemark® Implant.

Brånemark®	3.5 mm	4.1 mm	5.1 mm
Ti-Base K-Series	K 800	K 810	K 820
Abutment screw	K 60	K 61	K 62

The Abutment Ti-Base N-Series is compatible with Straumann SynOcta® Implant.

Straumann SynOcta®	3.5 mm	4.8 mm	6.5 mm
Ti-Base N-Series	N 800	N 810	N 820
Abutment screw	N 60	N 62	N 62

The Abutment Ti-Base R-Series is compatible with Zimmer Tapered Screw-vent® (Sulzer) Implant.

Zimmer (Sulzer) Tapered Screw-vent®	3.5 mm	4.5 mm	5.7 mm
Ti-Base R-Series	R 800	R 810	R 820
Abutment screw	R 60	R 60	R 60

The Abutment Ti-Base S-Series is compatible with Astra Tech OsseoSpeed® Implant.

Astra Tech Osseo-Speed®	3.5 / 4.0 mm	4.5 / 5.0 mm
Ti-Base S-Series	S 800 / 810	S 820
Abutment screw	S 60	S 61

The Abutment Ti-Base T-Series is compatible with Dentsply-Friadent Frialit® Implant.

Dentsply-Friadent Frialit®	3.4 mm	3.8 mm	4.5 mm	5.5 mm
Ti-Base T-Series	T 800	T 805	T 810	T 820
Abutment screw	T 60	T 60	T 60	T 60



The Abutment Ti-Base H-Series is compatible with Biomet 3i Osseotite® Certain®.

Biomet 3i Osseotite® Certain®	3.4 mm	4.1 mm	5.0 mm
Ti-Base H-Series	H800	H810	H 820
Abutment screw	H 60	H 60	H 60

The Abutment Ti-Base L-Series is compatible with Straumann® Bone Level Implant.

Straumann® BoneLevel	3.3 mm	4.1 / 4.8 mm
Ti-Base L-Series	L 800	L 810
Abutment screw	L 60	L 61

The Abutment Ti-Base F-Series is compatible with Nobel Biocare NobelActive™ Implant.

NobelActive™	3.5 mm	4.3 / 5.0 mm
Ti-Base F-Series	F 800	F 810
Abutment screw	F 60	F 61

Each Ti-Base is delivered with an abutment screw for fixation on the implant. The article number is the order number.

#### COMPOSITION:

Ti-Base and Abutment Screw: Ti6Al4V, medical grade 5, ASTM 136  
Scan Body: Polyether-ether-ketone, PEEK

#### CONTRAINDICATION:

The Ti-Bases of each Series can only be combined with the matching implant, e.g. the E-Series shall be combined exclusively with Replace select® implants. They cannot be combined with implants of a different implant type or manufacturer. The diameter of the Ti-base must correspond in size to the used implant in order to prevent a peri-implant tissue irritation.

The Ti-Bases are indicated for single use only. If they are used multiple times, they might damage the implants.

For fixation of the Ti-Bases on the implant, the correct torque force, recommended by the implant manufacturer, has to be considered carefully to avoid the damage of the implant-bone connection.

Mechanical treatment of the connection part of the Ti-Base will damage the correct fitting of the Ti-Base on the implant.

Ncm	Abutment				
20	H-Series				
25	S-Series	T-Series	F-Series		
30	R-Series				
35	I-Series	K-Series	N-Series	E-Series	L-Series

Mechanical treatment of the connection part of the Ti-Base will damage the correct fitting of the Ti-Base on the implant.

#### HANDLING METHOD FOR FURTHER PROCESSING:

The following instruction describes possible steps in the Laboratory for further processing to design the prosthetic components. The ceramic crown is an example for a possible aid prosthetic.

Ceramic crown: Milled with CAD/CAM-machines for zirconium oxide- or aluminum oxide- ceramic according to the anatomic form of a crown.

The ceramic crown shall be grinded or polished with diamond instruments and with minimal pressure and water-cooling. The minimal thickness shall be 0.5 mm. Sharpe edges must be avoided.

#### VENEERING:

Crowns shall be veneered with appropriate ceramics before cementing onto the Ti-Base. The instructions for use of the ceramic manufacturers have to be considered.

Before cementing the Ti-Base and the ceramic crown: Sandblasting the contact surfaces with Al<sub>2</sub>O<sub>3</sub>, 50 µm, 2 bar and ensure intensive cleaning of dust and grease.

It is recommended to protect the connection part of the Ti-Base with an implant analog during handling.

#### CEMENTING:

It is recommended to cement the ceramic abutment onto the Ti-Base with Panavia® F2.0 (Kuraray) with RelayUnicem® (3M Espe) or other equivalent cements. The instructions for use of the cements shall be followed carefully.

The Ti-Base shall be fixed onto an implant analog with the abutments screw. The head of the screw has to be covered with wax or resin. The mixed cement is applied onto the contact part of the Ti-Base. The abutment is pressed onto the Ti-Base. The final position is evaluated by slight rotation. The gap between abutment and the Ti-Base must be as small as possible. Remaining cement shall be removed immediately.

#### POLISHING:

After hardening the remaining cement shall be removed with rotating silicon instruments. The cement inside the screw channel has to be removed carefully.



**SCAN BODY INDICATIONS:**

The Scan Body can be used as auxiliary component, to determine the exact position and insertion angle of the implant. This is helpful for the further process to design the prosthetics.

**SCAN BODY:**

For the CAD/CAM scanning of the model, the Scan Body is used to indicate the position of the implant. The size of the Scan Body shall be corresponding to the original Implant system, implant diameter and Ti-Base Series. The chamfer of the Scan Body prevents the rotation of the ceramic abutment. The Scan Body is fixed on the implant analogue with the abutment screw. After correct positioning, there is no gap visible between implant and Scan Body. Rotation of the Scan Body is impossible.

**ARTICLE NO.:**

The article number of Scan Body and Titanium Base is a combination of the code for the Series: E, I, K, N, R, S, T, H, L and F (→ X) with the code W for Scan Body and for Titanium Base.

Titanium-Base	X 800	X 805	X 810	X 820	X 830
Scan-Body	X 00W	X 05W	X 10W	X 20W	X 30W

**WARNING:**

**Safety hint:** metal dust is harmful to your health. When milling and sandblasting use a suction extraction system and a breathing mask.

**SECONDARY EFFECTS:**

Allergies to the alloy or contents of the alloy or electrochemically based reactions may very rarely occur.

**REACTIONS:**

In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

**WARRANTY:**

10 Years on the mechanical stability of the Ti-Base, if it was processed according to the Instruction for use. Whether given verbally, in writing or by practical instructions, our recommendation for use is based upon own experience and trials and can only be considered as standard values. Our products are subject to a constant further development. Therefore alternations in construction and composition are reserved.

**CLEANING, DISINFECTION AND STERILIZATION:**

The nt-trading abutments and screws of the series E, I, K, N, R, S and T are supplied in non-sterile condition. The components should be cleaned, disinfected and in specific clinical procedures and cases be sterilized, prior to use after they are received from the dental laboratory (no liability on disregard). Effective cleaning and disinfection is an indispensable requirement for effective sterilization of the abutments. Prior to sterilization, please keep implants and screws clean when handling in the laboratory and operatory.

Additionally, please pay attention to legal regulations valid for your local areas as well as to the hygienic instructions of your dental practice. This applies particularly to the different guidelines regarding the inactivation of prions.

**1. Pre-disinfection** (avoidance of cross contaminations)

Place the abutments and screws in a germicidal bath\* immediately after use. Remove all residues and disassemble demountable products.

**2. Cleaning**

Please use distilled water and neutral cleaning agents\* only. The internal irrigation tube has to be cleaned with a Miller needle and must be rinsed with distilled water at the beginning and end of the exposure time using a disposable syringe (min. 10 ml). The products must be cleaned with a plastic instrument cleaning brush and then rinsed with distilled water. Please control all products after cleaning in order to avoid either damaging or corrosion. Damaged products must be replaced.

**3. Rinsing and Drying**

After removal of the products from the germicidal bath, all components must be rinsed 3 times with distilled water (e.g. Aqua purificata). Please dry all components thoroughly with a lint-free disposable cloth. For the cleaning of the internal irrigation tube oil-free compressed air is mandatory. Please re-check all parts for damage or corrosion afterwards.

**4. Disinfection**

We recommend a high level disinfectant such as, Cidex OPA (Johnson & Johnson) for disinfection of the abutments and screws.

- a) Soak the abutments and screws in the disinfectant solution for the required amount of time. See instructions for use of Cidex OPA.
- b) Remove the abutments and screws from the disinfectant solution.
- c) Rinse at least three times with highly purified water.
- d) Air dry and package the abutments and screws immediately.

**5. Sterilization:**

If no sterilization device is available in the laboratory, this information should be forwarded to the dentist so proper sterilization can occur. Please use only validated sterilization procedures for the sterilization of the abutments and screws. Other sterilization procedures must not be used.

**Reusability:**

You may only sterilize the abutments one time. In case of inadvertent contamination, you may re-sterilize one time after cleaning and disinfection.

**Steam sterilization:**

- fractionated vacuum procedure or gravity procedure (with sufficient product drying)
- steam sterilizer according to ISO 17665: 2006 or EN 13060 and EN 285 respectively or equivalent national standards
- validated according to EN ISO/ANSI AAMI 17665 (in past: EN 554/ANSI AAMI ISO 11134) (valid IQ/OQ) (commissioning and product specific performance qualification)
- sterilization time 20 minutes at 121 °C (250 °F) (listed exposure times are at sterilization temperature)

**7. Storage**

Store the sterilized parts dry and dust-free at room temperature.

\* Please observe all manufacturers guidelines for disinfection and cleaning agents with special regard to the concentration, exposure time and temperature. Only neutral disinfection solutions without chlorine, ammonia and aldehydes and with a proven effectiveness against HBV, HCV, HDV and HIV must be used. The products have to meet the respective national regulations for disinfectants. If disinfectants containing aldehydes are used, this might lead to a possible fixation of proteins. Please use only freshly prepared solutions.

Rev. A/ 2012-03-20

**MANUFACTURER:**

nt-trading GmbH & Co KG  
Essostrasse 16  
76187 Karlsruhe  
Germany  
Tel: +49 - 721 - 91 54 71 - 60  
Fax: +49 - 721 - 91 54 71 - 61  
E-mail: info@nt-trading.com

Products indicated with ® are registered brand names of the manufacturers.

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 5

---

**Appendix 5 – NT-Trading Ti-Base 510(k) Summary**

510(k) Summary

K111935

FEB 17 2012

Submitter Name: NT-Trading GmbH & Co. KG  
Submitter Address: Eссоstrasse 16  
76187 Karlsruhe  
Germany

Phone Number: +49-721-915471 60  
Fax Number: +49-721-915471 61

Contact Person: Dirk Jahn

Date Prepared: June 29, 2011

Device Trade Name: Ti-Base Abutment  
2-CONnect Abutment

Common Name: Dental Abutments

Classification Name, Number & Product Code: Abutment, Implant, Dental, Endosseous  
872.3630  
NHA

Predicate Devices: (K100152) Sirona Dental Systems Sirona Dental CAD/CAM System, (K083871) Atlantis™ Straumann Bone Level Abutment, (K093483) Atlantis™ Abutment for Nobel Active Implant, (K072642) Biomet 3I Dental Abutments And Restorative Components, (K990342) synOcta® Prosthetics, (K080239) P.004 Abutments, (K072570) NobelActive™ Multi Unit Abutment

Device Description and Statement of Intended Use

The Ti-Base Abutment is a premanufactured prosthetic component supplied in two parts, the abutment and screw, for fixation onto dedicated endosseous dental implants and is intended for use as an aid in prosthetic rehabilitation.

The 2-CONnect Abutment consists of 1 Abutment with screw (for fixation of abutment to the implant) and 1 titanium cap with 1 tiny screw (fixed into the hollow Abutment screw). The cap on top fits exactly to the abutment-geometry and does not have a rotation fixation, so it is easier to work with (not indicated for single crowns but strictly for bridges). The 2-CONnect is intended for use as an aid in prosthetic rehabilitation.

The NT-Trading Ti-Base and 2-CONnect is compatible with commercially available dental CAD/CAM systems, such as 3Shape, Exocad, Dental Wings. Such systems must be validated by the user.

Indication for use:  
**Ti-Base Abutments:** The devices covered by this submission are abutments which are placed into a dental implant to provide support for a prosthetic restoration.

The Ti-Base abutments are intended for use to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support single and multiple tooth prosthesis, in the mandible or maxilla. The prosthesis can be cement retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.

The Ti-Base abutments are indicated for use with the following implant systems:

---

Section 5.0: 510(k) Summary

- Nobel Biocare® Replace Select®
- Nobel Biocare NobelActive™
- Biomet 3i® Osseotite®
- Biomet 3i® Osseotite® Certain®
- Nobel Biocare Branemark®
- Straumann® synOcta®
- Straumann® Bone Level®
- Zimmer® Tapered Screw-vent®
- Astra Tech OsseoSpeed®
- Dentsply-Friadent® Frialit®

**2-CONnect Abutments:** 2-CONnect abutment is indicated for use to provide support for prosthetic restorations such as bars and bridges. The 2-CONnect abutments can be used in multiple tooth restorations. The 2-CONnect abutment can be used together with cemented bridges and bar constructions for functional and aesthetical reconstruction.

The 2-CONnect abutments are indicated for use with the following implant systems:

- Nobel Biocare® Replace Select®
- Straumann® synOcta®
- Straumann® BoneLevel®

Summary of Technological  
Characteristics

The proposed Ti-Base abutments and 2-CONnect abutments are substantially equivalent to the currently cleared devices. They are substantially equivalent in intended use, material and connection interfaces to the implants are identical for each individual diameter and connection type. Comparison Demonstrating Substantial Equivalence follows at the end of this section.

**Testing Summary**

In order to demonstrate compatibility of Ti-Base and 2-CONnect abutments to each implant system, fatigue testing was performed according to ISO 14801 Dentistry-Implants-Dynamic fatigue test for endosseous implants. Testing was performed on the abutments in this submission with the implants that they are intended to fit. See section 18.

Conclusion

The information discussed above demonstrates that the NT-Trading Ti-Base Dental Abutments and 2-CONnect Abutments are substantially equivalent to the predicate devices.

Declarations

- This summary includes only information that is also covered in the body of the 510(k).
- This summary does not contain any puffery or unsubstantiated labeling claims.
- This summary does not contain any raw data, i.e., contains only summary data.
- This summary does not contain any trade secret or confidential commercial information.
- This summary does not contain any patient identification information.

---

Section 5.0: 510(k) Summary

Summary of Technical Characteristics

Feature	Ti-Base and 2-CONNECT	Sirona Dental Systems Sirona Dental CAD/CAM System	Atlantis™ Straumann Bone Level Abutment	Atlantis™ Abutment for Nobel Active Implant	Biomet 3i Dental Restorations And Restorative Components	synOxide® Prosthetics	P.004 Abutments	Nobel/Active™ Multi Unit Abutment
510(k) Number		K100152	K083871	K093483	K072642	K990342	K080239	K072570
Manufacturer	Ni-Trading GmbH & Co. KG	Sirona Dental Systems GmbH	Astra Tech Inc.	Astra Tech Inc.	Biomet 3i, Inc.	Straumann® USA	Straumann® Manufacturing, Inc	Nobel Biocare® AB
Classification # & Product Code	872.3630 NHA	872.3630 NHA	872.3630 NHA	872.3630 NHA	872.3630 NHA	872.3630 NHA	872.3630 NHA	872.3630 NHA
Intended Use	<p>Ti-Base Abutments: The devices covered by this submission are placed into a dental implant to provide support for a prosthetic restoration.</p> <p>The Ti-Base abutments are intended for use to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support single and multiple tooth prostheses, in the mandible or maxilla. The prosthesis can be cement retained to the abutment. The abutment screw is</p>	<p>The Sirona Dental CAD/CAM System is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations. The system consists of three major parts: TiBase, InCoris mesostructure, and CAD/CAM software. Specifically, the InCoris mesostructure and TiBase components make up a two-piece abutment which is used in conjunction with endosseous dental implants to restore the function</p>	<p>The devices covered by this submission are placed into a dental implant to provide support for a prosthetic reconstruction. The Atlantis Abutment is intended for use with an endosseous implant to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support single and multiple tooth prostheses, in the mandible or maxilla. The prosthesis can be cement retained to the abutment. The abutment screw is intended to the endosseous implant.</p>	<p>The Atlantis Abutment is intended for use with an endosseous implant to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support single and multiple tooth prostheses, in the mandible or maxilla. The prosthesis can be cement retained to the abutment. The abutment screw is intended to the endosseous implant.</p>	<p>Dental Restorations and Overdenture Bars are intended for use as an accessory to endosseous dental implants to support a prosthetic device in a partially or edentulous patient. These are intended for use to support single and multiple tooth prostheses, in the mandible or maxilla. The prostheses can be cement retained to the abutment.</p>	<p>ITI Dental implants are intended to be placed in the maxillary and/or mandibular arch to support crowns, bridges or overdentures in edentulous or partially edentulous patients. The prosthetic accessories to dental implants are used either in the process of fabricating the prosthetic restoration for the implant or as part of the prosthetic restoration.</p>	<p>Abutments are placed into dental implants to provide support for prosthetic restorations such as crowns, bridges and overdentures. Abutments can be used in single tooth replacements and multiple tooth restorations. The subject abutments are for permanent screw-retained bridges and bar-retained implant-borne dentures. Permanent copings are intended to</p>	<p>Nobel Biocare's Multi-Unit is a premanufactured prosthetic component directly connected to endosseous dental implants and is intended for use as an aid in prosthetic rehabilitation.</p>

<p>intended to secure the abutment to the endosseous implant.</p> <p>The Ti-Base abutments are indicated for use with the following implant systems:</p> <ul style="list-style-type: none"> <li>• Nobel Biocare® Replace Select®</li> <li>• Nobel Biocare NobelActive™</li> <li>• Biomet 3i®</li> <li>• Osseofite®</li> <li>• Biomet 3i®</li> <li>• Osseofite®</li> <li>• Certain®</li> <li>• Nobel Biocare Branemark®</li> <li>• Straumann® synOcta®</li> <li>• Straumann® Bone Level®</li> <li>• Zimmer® Tapered Screw-vent®</li> <li>• Astra Tech OsseoSpeed®</li> <li>• Dentsply-Friadent® Frialfit®</li> </ul> <p><b>2-CONNECT Abutments:</b>                  2-CONNECT abutment is indicated for use to provide support for prosthetic restorations such as bars and bridges. The 2-CONNECT abutments can be</p>	<p>and aesthetics in the oral cavity. The InCoris mesostructure may also be used in conjunction with the Camlog Titanium base CAD/CAM (types K2244, xxxx) (K083496) in the Camlog Implant System. The CAD/CAM software is intended to design and fabricate the InCoris mesostructure. The InCoris mesostructure and TiBase two-piece abutment is compatible with the following implants systems:</p> <ul style="list-style-type: none"> <li>• Nobel Biocare Replace</li> <li>• Nobel Biocare Branemark</li> <li>• Friadent Xive</li> <li>• Biomet 3i</li> <li>• Osseofite</li> <li>• Astra Tech Osseospeed</li> <li>• Zimmer Tapered Screw-Vent</li> <li>• Straumann SynOcta</li> </ul>	<p>retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.</p>	<p>serve as a base for multi-unit bar or bridge restorations. Temporary Copings are intended to serve as a base for temporary restorations for up to 6 month. Protective Caps are intended to protect the outer configuration of the abutment and to maintain and condition the contours of the soft tissue during the healing phase for up to 6 months.</p>
---	---	---	--

	used in multiple tooth restorations. The 2-CONNECT abutment can be used together with cemented bridges and bar constructions for functional and aesthetical reconstruction. The 2-CONNECT abutments are indicated for use with the following implant systems:								
	<ul style="list-style-type: none"> <li>• Nobel Biocare®</li> <li>• Replace Select®</li> <li>• Straumann®</li> <li>• synOcta®</li> <li>• Straumann®</li> <li>• BoneLevel®</li> </ul>	Same	Same	Same	Same	Same	Same	Same	Same
Abutment Diameter min.	3.5 mm	Same	Same	Same	Same	Same	Same	Same	Same
Abutment Diameter max.	6.5 mm	Same	Same	Same	Same	Same	Same	Same	Same
Abutment Height	Ti-Base: 4 mm 2-CONNECT: 2.3 / 4.3 mm	Same	4 / 5.5 mm	6.6 mm	7.0	1.5 / 6.0	1.5 / 6.0	1.5 / 6.0	1.0 / 5.5 mm
Mode of Action	Screw-retained or cement retained	Screw-retained or cement retained	Screw-retained or cement retained	Screw-retained or cement retained	Screw-retained or cement retained	Screw-retained or cement retained	Screw-retained or cement retained	Screw-retained	Screw retained
Reusable	No	No	No	No	No	No	No	No	No
Material	Ti-6Al-4V	Ti-6Al-4V	Ti-6Al-4V ELI	Ti-6Al-4V ELI	Ti-6Al-4V ELI	Ti-6Al-4V ELI	Ti-6Al-4V	Titanium, Titanium alloy	Titanium Alloy

For the reasons stated above, we believe a determination of substantial equivalence between the Ti-Base and 2-CONNECT and these predicate devices is appropriate.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room –WO66-G609  
Silver Spring, MD 20993-0002

NT-Trading GmbH & Company AG  
C/O Mr. William Greenrose  
President  
Qserve America, Inc.  
220 River Road  
Claremont, New Hampshire 03743

FEB 17 2012

Re: K111935  
Trade/Device Name: Ti-Base for Individual milled Zirconium Abutment, 2-CONnect  
Abutment for Bridges and Bars  
Regulation Number: 21 CFR 872.3630  
Regulation Name: Endosseous Dental Implant Abutment  
Regulatory Class: II  
Product Code: NHA  
Dated: November 28, 2011  
Received: February 14, 2012

Dear Mr. Greenrose:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.



Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Anthony D. Watson, B.S., M.S., M.B.A.  
Director  
Division of Anesthesiology, General Hospital,  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K111935

Device Name: Ti-Base for individual milled Zirconium Abutment, 2-CONnect Abutment for Bridges and Bars

Indications For Use:

**Ti-Base for individual Zirconium Abutments:** The devices covered by this submission are abutments which are placed into a dental implant to provide support for a prosthetic restoration.

The Ti-Base for individual Zirconium Abutments are intended for use to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support single and multiple tooth prosthesis, in the mandible or maxilla. The prosthesis can be cement retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.

The Ti-Base abutments are indicated for use with the following implant systems:

- Nobel Biocare® Replace Select®
- Nobel Biocare NobelActive™
- Biomet 3i® Osseotite®
- Biomet 3i® Osseotite® Certain®
- Nobel Biocare Branemark®
- Straumann® synOcta®
- Straumann® Bone Level®
- Zimmer® Tapered Screw-vent®
- Astra Tech OsseoSpeed®
- Dentsply-Friadent® Frialit®

**2-CONnect Abutment for Bridges and Bars:** 2-CONnect Abutment for Bridges and Bars is indicated for use to provide support for prosthetic restorations such as bars and bridges. The 2-CONnect abutments can be used in multiple tooth restorations. The 2-CONnect abutment can be used together with cemented bridges and bar constructions for functional and aesthetical reconstruction.

The 2-CONnect abutments are indicated for use with the following implant systems:

- Nobel Biocare® Replace Select®
- Straumann® synOcta®
- Straumann® BoneLevel®

Prescription Use  X   
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
(Division Sign-Off)

Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K111935



Attachment 3  
02.11.2002

510(k) SUMMARY

K014174

A. Submitter's Name and Address

Barnstead/Thermolyne Corp.  
P.O. Box 797  
Dubuque, IA 52004

FEB 22 2002

B. Contact Person

Mia M. Ware  
Regulatory Affairs Specialist  
563-556-2241 Ext. 485  
Fax: 563-557-0612

C. Establishment Registration Number of Submitter

1950043

D. Contract Manufacturing Facility

Not Applicable

E. Device Name

Proprietary Name: Harvey® PV Dry  
Common Name: Steam Sterilizer  
Classification Name: Steam Sterilizer

F. Device Classification

Class II §880.6880

G. Action Taken to Comply with Section 514 of the Act

The Agency has recognized the ANSI/AAMI ST55, Tabletop Steam Sterilizers, and the FDA Guidance on Premarket Notification Submissions for Sterilizers Intended for Use in Health Care Facilities, March 1993, and its addendum, dated September, 1995. Conformance or variance with these standards is described on the following pages.

H. Reason for 510(k) Submission

- Initial Product Introduction
- New Model for Product-line Extension
- Initial Import into the USA
- Other (Include in Part IV an explanation referenced to Part I. H.)

I. Predicate Device: Harvey® MC10 Steam Sterilizer, K924955

Attachment 3  
02.11.2002

Description of the Device: The Harvey® PV Dry is a pre/post-vacuum table top steam sterilizer. Its enclosure is made of painted metal and door cover is made of plastic. The outside dimensions are 15.5"Hx19"Wx23"D, it weighs 115lbs and the chamber is 10" in diameter. The chamber is made from stainless steel rated at 45 psig to comply with the American Society of Mechanical Engineers (ASME) Pressure Vessel Code.

Intended Use of the Device: The Harvey® PV Dry is a pre/post-vacuum table top steam sterilizer that is designed for use in medical and dental offices, hospitals, clinics, and other facilities where a variety of materials require sterilization. Intended for sterilization of wrapped or unwrapped instruments, dental handpieces, and linen packs.

Technological Characteristics: The Harvey® PV Dry has very similar technological characteristics as the predicate device, Harvey® MC10 Steam Sterilizer. The Harvey® PV Dry is a table-top autoclave steam sterilizer with pre-vacuum air removal and sterile post-vacuum drying, similar to larger hospital type units.

The Harvey® PV Dry provides superior performance for dental handpieces and difficult to penetrate fabric packs and superior and rapid drying for all loads. With sterile vacuum drying, the Harvey® PV Dry is recommended for medical and dental offices where packs, bagged instruments, and wrapped instrument sets, such as surgical kits are stored for later usage. The integral sterile drying in a 30-40 minute cycle, gives assurance of sterile conditions of the instruments at the time of delivery to the operatory. Sterile vacuum drying eliminates the need for extended drying cycles using non-sterile air.

Some features unique to the Harvey® PV Dry also demonstrate safety as well as the efficiency of the unit. The interlocking device on the door prevents manual opening until the chamber pressure decreases to near atmospheric pressure. Other technologies which help in the convenience and ease of operation of the Harvey® PV Dry include indicator lights to tell the user when the waste tank needs to be emptied, and controls preventing the running of a cycle until it is drained, indicator lights telling the user when the water supply needs to be replenished, and cycle parameter display. In the predicate device comparison matrix in Section II, there are detailed differences and similarities between the Harvey® PV Dry and its predicate device, Harvey® MC10 Steam Sterilizer.

Non-clinical Testing: Validation studies were conducted by SPS Medical Supply Corporation located at 6789 West Henrietta Road, Rush, NY 14543. Validation testing is in accordance with ANSI/AAMI ST55, ANSI/AAMI ST37. Since neither standard specifically addresses handpiece cycle validation, additional handpiece cycle validation testing was performed with the protocol being written in accordance with FDA Draft Guidance Document on Dental Handpieces, 1995, and the Guidance on Premarket Notification Submissions for Sterilizers Intended for Use in Health Care Facilities, 1993. Successful sterilization was achieved in all three validations and the declarations of conformity to the consensus standards are located in Section III. Testing results and raw data are contained in the product's device master record.

Conclusion: It is Barnstead/Thermolyne's conclusion that the Harvey® PV Dry is substantially equivalent to its predicate device, the Harvey® MC10 Steam Sterilizer. Based on the validation results and information submitted, the Harvey® PV Dry provides effective sterilization of the indicated medical and dental materials.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

FEB 22 2002

Barnstead/Thermolyne Corporation  
C/O Mr. Reiner Krumme  
Responsible Third Party  
TUV Rheinland of North America, Inc  
12 Commerce Road  
Newton, Connecticut 06470

Re: K014174

Trade/Device Name: Harvey PV Dry Sterilizer, Model # ST127325; ST127320

Regulation Number: 880.6880

Regulation Name: Steam Sterilizer

Regulatory Class: II

Product Code: FLE

Dated: February 12, 2002

Received: February 13, 2002

Dear Ms. Krumme:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements


Page 2 – Mr. Krumme

of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4618. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



Timothy A. Ulatowski  
Director  
Division of Dental, Infection Control  
and General Hospital Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

Attachment #  
02.11.2002

### INDICATIONS FOR USE STATEMENT

510(k) Number : 014174

Device Name: Harvey® PV Dry

**Indications for Use:**

The Harvey® PV Dry is a prevacuum and post vacuum drying sterilizer intended to provide sterilization of medical and dental instruments. It is intended to provide sterilization of wrapped or unwrapped metal instruments, surgical devices and other heat stable devices in pouches, surgical packs, dental handpieces, and linen packs. Deionized or distilled water is required for operation of the Harvey® PV Dry.


The sterilizer has four standard cycles: Unwrapped Instruments (135C for 3 minutes); Wrapped Instruments (135C for five minutes); Packs (121C for 30 minute); Special, for dental handpieces (5 minutes at 134C). There is also an accessible Bowie-Dick test cycle for routine testing of the steam penetration capability into packs. The parameters are fixed for each cycle.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Optional Format 3-10-98)



(Division Sign-Off)

Division of Dental, Infection Control,

General Hospital Devices

Device number \_\_\_\_\_

1014174



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 7

---

**Appendix 7 – (b)(4) TS/CCI Test Report (b)(4) TS/CCI**





















































































**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**


Appendix 10

---

(b)(4) Trade Secret Process - Product Specs









**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 11

---

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 11

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 12

---

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 12

---

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Appendix 13

---

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 13

(b)(4) Trade Secret Process - Product Specs





**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Appendix 14

---

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 14

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Appendix 15

---

(b)(4) Trade Secret Process - Product Specs























January 23, 2014

U.S. Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center – W066-G609  
10903 New Hampshire Avenue  
Silver Spring, Maryland 20993-0002

**Subject: Response to Telephone Holds for K132219 Dated November 14, 2013 and  
November 19, 2013**

Dear Sir or Madam:

(b)(4) Trade Secret Process - Product Specs



(b)(4) Trade Secret Process - Product Specs




(b)(4) Trade Secret Process - Product Specs

A large, solid black rectangular redaction box covers the majority of the page content, starting below the header and ending above the footer. The text '(b)(4) Trade Secret Process - Product Specs' is written in red at the top left corner of this redacted area.

(b)(4) Trade Secret Process - Product Specs

A large, solid black rectangular redaction box covers the majority of the page's content, starting below the header and ending above the footer. The text '(b)(4) Trade Secret Process - Product Specs' is printed in red at the top left corner of this redacted area.





(b)(4) Trade Secret Process - Product Specs



Sincerely,



Jennifer M. Jackson, MS

Senior Regulatory Affairs Project Manager

Enclosures

---



# **Traditional 510(k) Submission**

**Straumann® Variobase™ Abutments**

CDRH Premarket Review Submission Cover Sheet

---

## **2 CDRH Premarket Review Submission Cover Sheet**

The CDRH Premarket Review Submission Cover Sheet begins on the next page.

**CDRH PREMARKET REVIEW SUBMISSION COVER SHEET**

Date of Submission 01-23-2014	User Fee Payment ID Number <b>(b)(4) Trade</b>	FDA Submission Document Number (if known) K132219
----------------------------------	---	--

**SECTION A TYPE OF SUBMISSION**

<b>PMA</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Premarket Report <input type="checkbox"/> Modular Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment <input type="checkbox"/> Licensing Agreement	<b>PMA &amp; HDE Supplement</b> <input type="checkbox"/> Regular (180 day) <input type="checkbox"/> Special <input type="checkbox"/> Panel Track (PMA Only) <input type="checkbox"/> 30-day Supplement <input type="checkbox"/> 30-day Notice <input type="checkbox"/> 135-day Supplement <input type="checkbox"/> Real-time Review <input type="checkbox"/> Amendment to PMA & HDE Supplement <input type="checkbox"/> Other	<b>PDP</b> <input type="checkbox"/> Original PDP <input type="checkbox"/> Notice of Completion <input type="checkbox"/> Amendment to PDP	<b>510(k)</b> <input checked="" type="checkbox"/> Original Submission: <input checked="" type="checkbox"/> Traditional <input type="checkbox"/> Special <input type="checkbox"/> Abbreviated (Complete section I, Page 5) <input type="checkbox"/> Additional Information <input type="checkbox"/> Third Party	<b>Request for Feedback</b> <input type="checkbox"/> Pre-Submission <input type="checkbox"/> Informational Meeting <input type="checkbox"/> Submission Issue Meeting <input type="checkbox"/> Day 100 Meeting <input type="checkbox"/> Agreement Meeting <input type="checkbox"/> Determination Meeting <input type="checkbox"/> Study Risk Determination <input type="checkbox"/> Other (specify):
<b>IDE</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement	<b>Humanitarian Device Exemption (HDE)</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment	<b>Class II Exemption Petition</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information	<b>Evaluation of Automatic Class III Designation (De Novo)</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information	<b>Other Submission</b> <input type="checkbox"/> 513(g) <input type="checkbox"/> Other (describe submission):

Have you used or cited Standards in your submission?  Yes  No (If Yes, please complete Section I, Page 5)

**SECTION B SUBMITTER, APPLICANT OR SPONSOR**

Company / Institution Name Straumann USA, LLC	Establishment Registration Number (if known) 1222315		
Division Name (if applicable)	Phone Number (including area code) 978-747-2509		
Street Address 60 Minuteman Road	FAX Number (including area code) 978-747-0023		
City Andover	State / Province MA	ZIP/Postal Code 01810	Country USA
Contact Name Jennifer M. Jackson, MS			
Contact Title Senior Regulatory Affairs Project Manager		Contact E-mail Address jennifer.jackson@straumann.com	

**SECTION C APPLICATION CORRESPONDENT (e.g., consultant, if different from above)**

Company / Institution Name			
Division Name (if applicable)		Phone Number (including area code)	
Street Address		FAX Number (including area code)	
City	State / Province	ZIP Code	Country
Contact Name			
Contact Title		Contact E-mail Address	

**SECTION D1 REASON FOR APPLICATION - PMA, PDP, OR HDE**

<input type="checkbox"/> New Device <input type="checkbox"/> Withdrawal <input type="checkbox"/> Additional or Expanded Indications <input type="checkbox"/> Request for Extension <input type="checkbox"/> Post-approval Study Protocol <input type="checkbox"/> Request for Applicant Hold <input type="checkbox"/> Request for Removal of Applicant Hold <input type="checkbox"/> Request to Remove or Add Manufacturing Site	<input type="checkbox"/> Change in design, component, or specification: <input type="checkbox"/> Software / Hardware <input type="checkbox"/> Color Additive <input type="checkbox"/> Material <input type="checkbox"/> Specifications <input type="checkbox"/> Other ( <i>specify below</i> )	<input type="checkbox"/> Location change: <input type="checkbox"/> Manufacturer <input type="checkbox"/> Sterilizer <input type="checkbox"/> Packager
<input type="checkbox"/> Process change: <input type="checkbox"/> Manufacturing <input type="checkbox"/> Packaging <input type="checkbox"/> Sterilization <input type="checkbox"/> Other ( <i>specify below</i> )	<input type="checkbox"/> Labeling change: <input type="checkbox"/> Indications <input type="checkbox"/> Instructions <input type="checkbox"/> Performance Characteristics <input type="checkbox"/> Shelf Life <input type="checkbox"/> Trade Name <input type="checkbox"/> Other ( <i>specify below</i> )	<input type="checkbox"/> Report Submission: <input type="checkbox"/> Annual or Periodic <input type="checkbox"/> Post-approval Study <input type="checkbox"/> Adverse Reaction <input type="checkbox"/> Device Defect <input type="checkbox"/> Amendment
<input type="checkbox"/> Response to FDA correspondence:		<input type="checkbox"/> Change in Ownership <input type="checkbox"/> Change in Correspondent <input type="checkbox"/> Change of Applicant Address

Other Reason (*specify*):

**SECTION D2 REASON FOR APPLICATION - IDE**

<input type="checkbox"/> New Device <input type="checkbox"/> New Indication <input type="checkbox"/> Addition of Institution <input type="checkbox"/> Expansion / Extension of Study <input type="checkbox"/> IRB Certification <input type="checkbox"/> Termination of Study <input type="checkbox"/> Withdrawal of Application <input type="checkbox"/> Unanticipated Adverse Effect <input type="checkbox"/> Notification of Emergency Use <input type="checkbox"/> Compassionate Use Request <input type="checkbox"/> Treatment IDE <input type="checkbox"/> Continued Access	<input type="checkbox"/> Change in: <input type="checkbox"/> Correspondent / Applicant <input type="checkbox"/> Design / Device <input type="checkbox"/> Informed Consent <input type="checkbox"/> Manufacturer <input type="checkbox"/> Manufacturing Process <input type="checkbox"/> Protocol - Feasibility <input type="checkbox"/> Protocol - Other <input type="checkbox"/> Sponsor	<input type="checkbox"/> Response to FDA Letter Concerning: <input type="checkbox"/> Conditional Approval <input type="checkbox"/> Deemed Approved <input type="checkbox"/> Deficient Final Report <input type="checkbox"/> Deficient Progress Report <input type="checkbox"/> Deficient Investigator Report <input type="checkbox"/> Disapproval <input type="checkbox"/> Request Extension of Time to Respond to FDA <input type="checkbox"/> Request Meeting <input type="checkbox"/> Request Hearing
	<input type="checkbox"/> Report submission: <input type="checkbox"/> Current Investigator <input type="checkbox"/> Annual Progress Report <input type="checkbox"/> Site Waiver Report <input type="checkbox"/> Final	

Other Reason (*specify*):

**SECTION D3 REASON FOR SUBMISSION - 510(k)**

<input checked="" type="checkbox"/> New Device	<input type="checkbox"/> Additional or Expanded Indications	<input type="checkbox"/> Change in Technology
--	---	---

Other Reason (*specify*):

**SECTION E ADDITIONAL INFORMATION ON 510(K) SUBMISSIONS**

Product codes of devices to which substantial equivalence is claimed								Summary of, or statement concerning, safety and effectiveness information	
1	NHA	2		3		4		<input checked="" type="checkbox"/> 510 (k) summary attached <input type="checkbox"/> 510 (k) statement	
5		6		7		8			

Information on devices to which substantial equivalence is claimed (if known)

	510(k) Number	Trade or Proprietary or Model Name	Manufacturer
1	K120822	Straumann CARES Variobase Abutment NNC, RN, WN, NC, RC	Institut Straumann AG
2	K111935	Ti-Base Abutment	NT-Trading GmbH & Co. KG
3	K072055	Lava Frame, Lava Frame Shade	3M ESPE AG
4	K072071	P.004 RC Cementable Abutments	Institut Straumann AG
5	K111421	Sirona Dental CAD/CAM-System	Sirona Dental Systems GmbH
6	K120053	IPS e.max Press - Abutment Solutions	Ivoclar Vivadent AG

**SECTION F PRODUCT INFORMATION - APPLICATION TO ALL APPLICATIONS**

Common or usual name or classification name

Endosseous dental implant abutment

(b)(4) Trade Secret Process - Testing

FDA document numbers of all prior related submissions (regardless of outcome)

1	None	2		3		4		5		6	
7		8		9		10		11		12	

Data Included in Submission

Laboratory Testing       Animal Trials       Human Trials

**SECTION G PRODUCT CLASSIFICATION - APPLICATION TO ALL APPLICATIONS**

Product Code NHA	C.F.R. Section (if applicable) 21 CFR 872.3630	Device Class <input type="checkbox"/> Class I <input checked="" type="checkbox"/> Class II <input type="checkbox"/> Class III <input type="checkbox"/> Unclassified
Classification Panel Dental		

Indications (from labeling)

The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth restorations.

**SECTION E ADDITIONAL INFORMATION ON 510(K) SUBMISSIONS**

Product codes of devices to which substantial equivalence is claimed								Summary of, or statement concerning, safety and effectiveness information  <input type="checkbox"/> 510 (k) summary attached <input type="checkbox"/> 510 (k) statement
1		2		3		4		
5		6		7		8		

Information on devices to which substantial equivalence is claimed (if known)

	510(k) Number	Trade or Proprietary or Model Name	Manufacturer
1	K132209	IPS e.max CAD Abutment Solutions	Ivoclar Vivadent AG
2			
3			
4			
5			
6			

**SECTION F PRODUCT INFORMATION - APPLICATION TO ALL APPLICATIONS**

Common or usual name or classification name

	Trade or Proprietary or Model Name for This Device	Model Number
1		1
2		2
3		3
4		4
5		5

FDA document numbers of all prior related submissions (regardless of outcome)

1	2	3	4	5	6
7	8	9	10	11	12

Data Included in Submission  
 Laboratory Testing       Animal Trials       Human Trials

**SECTION G PRODUCT CLASSIFICATION - APPLICATION TO ALL APPLICATIONS**

Product Code	C.F.R. Section (if applicable)	Device Class
		<input type="checkbox"/> Class I <input type="checkbox"/> Class II <input type="checkbox"/> Class III <input type="checkbox"/> Unclassified
Classification Panel		

Indications (from labeling)

**Note:** Submission of the information entered in Section H does not affect the need to submit device establishment registration.

FDA Document Number (if known)

**SECTION H MANUFACTURING / PACKAGING / STERILIZATION SITES RELATING TO A SUBMISSION**

<input checked="" type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete		Facility Establishment Identifier (FEI) Number		<input checked="" type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Contract Manufacturer <input type="checkbox"/> Repackager / Relabeler	
Company / Institution Name Institut Straumann AG		Establishment Registration Number 9613348			
Division Name (if applicable)		Phone Number (including area code) 978-747-2509			
Street Address Peter Merian-Weg 12		FAX Number (including area code) 978-747-0023			
City Basel		State / Province		ZIP Code CH-4052	Country Switzerland
Contact Name Jennifer M. Jackson, MS		Contact Title Senior Regulatory Affairs Project Manager		Contact E-mail Address jennifer.jackson@straumann.com	

<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete		Facility Establishment Identifier (FEI) Number		<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Contract Manufacturer <input type="checkbox"/> Repackager / Relabeler	
Company / Institution Name		Establishment Registration Number			
Division Name (if applicable)		Phone Number (including area code)			
Street Address		FAX Number (including area code)			
City		State / Province		ZIP Code	Country
Contact Name		Contact Title		Contact E-mail Address	

<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete		Facility Establishment Identifier (FEI) Number		<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Contract Manufacturer <input type="checkbox"/> Repackager / Relabeler	
Company / Institution Name		Establishment Registration Number			
Division Name (if applicable)		Phone Number (including area code)			
Street Address		FAX Number (including area code)			
City		State / Province		ZIP Code	Country
Contact Name		Contact Title		Contact E-mail Address	




**SECTION I**

**UTILIZATION OF STANDARDS**

**Note:** Complete this section if your application or submission cites standards or includes a "Declaration of Conformity to a Recognized Standard" statement.

(b)(4) Trade Secret Process - Testing



**Please include any additional standards to be cited on a separate page.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 0.5 hour per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
1350 Piccard Drive, Room 400  
Rockville, MD 20850

*An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number.*

**SECTION I**

**UTILIZATION OF STANDARDS**

**Note:** Complete this section if your application or submission cites standards or includes a "Declaration of Conformity to a Recognized Standard" statement.

	Standards No.	Standards Organization	Standards Title	Version	Date
1	ISO 17664:2004	ISO	Sterilization of medical devices -- Information to be provided by the manufacturer for the processing of reesterilizable medical devices	2004	10-12-2012
2	ISO 14801:2007	ISO	Dentistry -- Implants -- Dynamic fatigue test for endosseous dental implants	2007	04-25-2012
3					
4					
5					
6					
7					

**Please include any additional standards to be cited on a separate page.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 0.5 hour per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
 Food and Drug Administration  
 Office of Chief Information Officer  
 Paperwork Reduction Act (PRA) Staff  
 1350 Piccard Drive, Room 400  
 Rockville, MD 20850

*An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number.*

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Indications for Use Statement

---

**4 Indications for Use Statement**

The Indications for Use Statement associated with this 510(k) is located on the following page in the required format.

## Indications for Use

510(k) Number (if known):

Device Name: Straumann® Variobase™ Abutments

Indications for Use:

The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth restorations.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF  
NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
510(k) Summary

---

## **5 510(k) Summary**

### **5.1 Submitter's Contact Information**

Straumann USA, LLC (on behalf of Institut Straumann AG)

60 Minuteman Road

Andover, MA 01810

Phone Number: 1-978-747-2509

Fax Number: 1-978-747-0023

Contact Person: Jennifer M. Jackson, MS

Date of Submission: 23-Jan-2014

### **5.2 Name of the Device**

Trade Name: Straumann® Variobase™ Abutments

Common Name: Dental Implant Abutment

Classification Name: Abutment, Implant, Dental, Endosseous

Regulation Number: §872.3630

### **5.3 Predicate Device(s)**

- K120822, Straumann CARES Variobase Abutment NNC, RN, WN, NC, RC (Institut Straumann AG)
- K111935, Ti-Base Abutment (NT-Trading GmbH & Co. KG)
- K072055, Lava™ Frame, Lava™ Frame Shade (3M ESPE AG)
- K072071, P.004 RC Cementable Abutments (Institut Straumann AG)
- K111421, Sirona Dental CAD/CAM-System (Sirona Dental Systems GmbH)
- K120053, IPS e.max® Press – Abutment Solutions (Ivoclar Vivadent AG)
- K132209, IPS e.max® CAD Abutment Solutions (Ivoclar Vivadent AG)

### **5.4 Device Description**

The Straumann® Variobase™ Abutments are pre-manufactured (stock) abutments, sometimes referred to as “Ti-bases”. Straumann® Variobase™ Abutments are available

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### 510(k) Summary

---

to fit Straumann® dental implant platforms NNC (Narrow Neck CrossFit®), RN (Regular Neck), WN (Wide Neck), NC (Narrow CrossFit®), and RC (Regular CrossFit®). A dental laboratory technician would design the corresponding coping and/or crown (the second component of the Variobase two-piece abutment) and/or prosthetic restoration in the dental laboratory via their preferred workflow for pressing, casting, or milling using either a burnout coping or STL model for open CAD software. The coping and/or crown would be manufactured via traditional laboratory methods for pressing or casting, or via validated Straumann milling.

#### 5.5 Intended Use

Straumann® Variobase™ Abutments are pre-manufactured (stock) abutments intended to be placed onto Straumann dental implants to provide support for customized prosthetic restorations such as crowns. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth restorations.

#### 5.6 Technological Characteristics

Straumann® Variobase™ Abutments are two-piece abutments consisting of a pre-manufactured (stock) abutment made from a titanium-aluminum-niobium alloy and a coping and/or crown which is designed in the dental laboratory by a dental technician and manufactured via traditional in-lab methods of pressing or casting, or via validated Straumann milling.

The Ti-base components of the Straumann® Variobase™ Abutments are identical to the Ti-base components of the Straumann predicate (K120822). The Ti-base components are also equivalent to the Ti-base components identified in K111935, K072055, K120053, and K132209.

The materials which may be used to manufacture the coping/crown component of the Straumann Variobase Abutments are identical to the identified predicate devices and include:

- Casting: Type 4 metals (ISO 22674)
- Pressing: IPS e.max® Press Ceramic (K120053)
- Milling: Polycon® ae (temporary restorations – K120822)

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### 510(k) Summary

---

Zerion® (K120822)

IPS e.max® CAD Ceramic (K132209)

Lava™ Frame, Lava™ Frame Shade (K072055)

### 5.7 Performance Testing

The material used in the manufacture of Straumann® Variobase™ Abutments is a titanium-aluminum-niobium alloy which meets the requirements of ISO 5832-11. Bench testing was performed with Polycon® ae and Zerion® to evaluate the fatigue load limits of the proposed Straumann® Variobase™ Abutments. Dynamic fatigue tests were conducted in accordance with the FDA guidance document *“Guidance for Industry and FDA Staff – Class II Special Controls Guidance Document: Root-form Endosseous Dental Implants and Endosseous Dental Abutments”*.

### 5.8 Conclusion

The documentation submitted in this premarket notification demonstrates that the Straumann® Variobase™ Abutments are substantially equivalent to the predicate devices and do not pose new issues of safety and effectiveness when used as labeled.

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Device Description

---

**10 Device Description**

(b)(4) Trade Secret Process - Testing



**10.1 Straumann® Variobase™ Abutments**

(b)(4) Trade Secret Process - Testing





**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Device Description

(b)(4) Trade Secret Process - Testing



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Device Description

---

Table 3 outlines the Straumann® Variobase™ Abutments and their corresponding article numbers:

(b)(4) Trade Secret Process - Testing



**10.2 Basal Screws**

(b)(4) Trade Secret Process - Testing



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Device Description

(b)(4) Trade Secret Process - Testing



**10.3 Accessories**

(b)(4) Trade Secret Process - Testing



# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

(b)(4) Trade Secret Process - Testing



### 10.4 Procedure

(b)(4) Trade Secret Process - Testing



#### 10.4.1 Restoration Design and Manufacturing

(b)(4) Trade Secret Process - Testing



##### 10.4.1.1 In-lab Casting

(b)(4) Trade Secret Process - Testing



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Device Description

(b)(4) Trade Secret Process - Testing



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Device Description

---

(b)(4) Trade Secret Process - Testing



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Device Description

(b)(4) Trade Secret Process - Testing



# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Substantial Equivalence Discussion

---

## 11 Substantial Equivalence Discussion

Within the meaning of the Medical Device Amendments Act of 1976, the proposed change to the Indications for Use for the Straumann® Variobase™ Abutments in this 510(k) premarket notification are substantially equivalent to the medical devices currently in commercial distribution listed below:

- K120822, Straumann CARES Variobase Abutment NNC, RN, WN, NC, RC (Institut Straumann AG)
- K111935, Ti-Base Abutment (NT-Trading GmbH & Co. KG)
- K072055, Lava™ Frame, Lava™ Frame Shade (3M ESPE AG)
- K072071, P.004 RC Cementable Abutments (Institut Straumann AG)
- K111421, Sirona Dental CAD/CAM-System (Sirona Dental Systems GmbH)
- K120053, IPS e.max® Press – Abutment Solutions (Ivoclar Vivadent AG)
- K132209, IPS e.max® CAD Abutment Solutions (Ivoclar Vivadent AG)

Straumann® Variobase™ Abutments are pre-manufactured (stock) abutments intended to be placed onto Straumann dental implants to provide support for customized prosthetic reconstructions such as crowns. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth restorations.

The NT-Trading Ti-Base Abutment is a pre-manufactured abutment supplied in two parts, the abutment and screw, for fixation onto dedicated endosseous dental implants and is intended for use as an aid in prosthetic restoration. As with the subject device, the coping/restoration is provided by the dental laboratory. The 510(k) Summary for the NT-Trading Ti-Base Abutment (Appendix 5) states “The Ti-Base is compatible with commercially available dental CAD/CAM systems, such as 3Shape, Exocad, and Dental Wings. Such systems must be validated by the user”.

The NT-Trading catalog and package insert are included in this submission in Appendices 3 and 4, respectively. The products that were cleared in premarket notification K111935 are outlined in Table 7 (the 510(k) Summary is included in Appendix 5). Specifically, the L-Series and N-Series abutments are compatible with implants of the Straumann Dental Implant System as shown in Table 7.



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Substantial Equivalence Discussion

(b)(4) Trade Secret Process - Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Substantial Equivalence Discussion

(b)(4) Trade Secret Process - Product Specs



## Traditional 510(k) Submission

### Straumann® Variobase™ Abutments

#### Substantial Equivalence Discussion

The table below provides a comparison matrix of the proposed and predicate devices (K120822):

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
<b>K Number</b>	<b>Straumann® Variobase™ Abutments Subject Submission</b>	<b>Straumann® CARES® Variobase™ Abutments (K120822)</b>	
<b>Indications for Use</b>	<p>The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth restorations.</p>	<p>The Straumann® CARES® Variobase™ Abutment is a two-piece dental abutment consisting of the Straumann® Variobase™ Abutment and the Straumann® CARES® Variobase™ Coping which is intended to be placed onto Straumann dental implants to provide support for prosthetic reconstruction such as crown and bridges. Straumann® CARES® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth and bridge restorations.</p> <p>The Straumann® CARES® Variobase™ Coping polycon® ae in combination with the Straumann® CARES® Variobase™ Abutment is indicated for temporary (up to 180 days) dental restoration of a Straumann dental implant.</p>	<p>The indications for use are being modified to allow Straumann to market the Straumann® Variobase™ Abutment as a stand-alone component. The dental laboratory would then manufacture the respective coping and/or prosthetic restoration using a burnout coping or STL model for open CAD software.</p>
<b>Ti-base Material</b>	Titanium-Aluminum-Niobium alloy (Ti-6Al-7Nb)	Titanium-Aluminum-Niobium alloy (Ti-6Al-7Nb)	Identical
<b>Abutment Diameter</b>	(b)(4) Trade Secret Process - Product Specs		Identical
<b>Abutment Height</b>			Identical

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
 Substantial Equivalence Discussion

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
K Number	Straumann® Variobase™ Abutments Subject Submission	Straumann® CARES® Variobase™ Abutments (K120822)	
Bonding Surface Area	(b)(4) Trade Secret Process - Product Specs		Identical
Coping/ Crow Material	(b)(4) Trade Secret Process - Product Specs		
Design Workflow	Wax-up or Open CAD	CARES® Visual	Equivalent
Manufacturing Workflow	(b)(4) Trade Secret Process - Product Specs		
Mode of Action	Screw-retained or cement retained	Screw-retained or cement retained	Identical
Reusable	No	No	Identical

**Table 8 - Comparison Matrix: Proposed Device versus Predicate Devices (K120822)**

## Traditional 510(k) Submission

### Straumann® Variobase™ Abutments

#### Substantial Equivalence Discussion

The table below provides a comparison matrix of the proposed and predicate devices (K111935):

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
<b>K Number</b>	<b>Straumann® Variobase™ Abutments Subject Submission</b>	<b>Ti-Base Abutment (K111935)</b>	
<b>Indications for Use</b>	<p>The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth restorations.</p>	<p>Ti-Base for individual Zirconium Abutments: The devices covered by this submission are abutments which are placed into a dental implant to provide support for a prosthetic restoration.</p> <p>The Ti-Base for individual Zirconium Abutments are intended for use to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support single and multiple tooth prosthesis, in the mandible or maxilla. The prosthesis can be cement retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.</p> <p>The Ti-Base abutments are indicated for use with the following implant systems:</p> <ul style="list-style-type: none"> <li>• Nobel Biocare® Replace Select®</li> <li>• Nobel Biocare NobelActive™</li> <li>• Biomet 3i® Osseotite®</li> <li>• Biomet 3i® Osseotite® Certain®</li> <li>• Nobel Biocare Branemark®</li> <li>• Straumann® synOcta®</li> <li>• Straumann® Bone Level®</li> <li>• Zimmer® Tapered Screw-vent®</li> <li>• Astra Tech OsseoSpeed®</li> <li>• Dentsply-Friadent® Frialit®</li> </ul>	<p>Equivalent – Both the subject and predicate devices are designed to interface with the Straumann Bone Level or the Straumann Tissue Level implants.</p>

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
 Substantial Equivalence Discussion

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
K Number	Straumann® Variobase™ Abutments Subject Submission	Ti-Base Abutment (K111935)	
Ti-base Material	Titanium-Aluminum-Niobium alloy (Ti-6Al-7Nb)	Titanium-Aluminum-Vanadium alloy (Ti-6Al-4V)	Equivalent
Abutment Diameter	(b)(4) Trade Secret Process - Product Specs		Equivalent
Abutment Height			Equivalent
Bonding Surface Area	(b)(4) Trade Secret Process - Product Specs		Unknown
Coping/ Cro Material			Not specified
Design Workflow	(b)(4) Trade Secret Process - Product Specs		

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
 Substantial Equivalence Discussion

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
<b>K Number</b>	<b>Straumann® Variobase™ Abutments Subject Submission</b>	<b>Ti-Base Abutment (K111935)</b>	
<b>Manufacturing Workflow</b>	<b>(b)(4) Trade Secret Process - Product Specs</b>		
<b>Mode of Action</b>	Screw-retained or cement retained	Screw-retained or cement retained	Identical
<b>Reusable</b>	No	No	Identical

**Table 9 - Comparison Matrix: Proposed Device versus Predicate Devices (K111935)**

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Substantial Equivalence Discussion

The table below provides a comparison matrix of the proposed and predicate devices (K120053/K072071 and K132209/K111421):

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
<b>K Number</b>	<b>Straumann® Variobase™ Abutments Subject Submission</b>	<b>IPS e.max® Press – Abutment Solutions (K120053) including P.004 RC Cementable Abutments (K072071)</b>	<b>IPS e.max CAD Abutment Solutions (K132209) including Sirona Dental CAD/CAM System (K111421)</b>	
<b>Indications for Use</b>	The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth restorations.	IPS e.max® Press Abutment Solutions is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations. IPS e.max Press Abutment Solutions is recommended for the fabrication of: - Hybrid abutments for single-tooth restorations - Hybrid abutment crowns for restorations The following Ti bases are intended to be used with IPS e.max Press Abutment Solutions: Straumann® Bone Level RC Ø4.1 mm or Ø4.8 mm (K062129) RC Cementable Abutment D 5.0 – 6.5 mm, GH 1.0 – 3.0 mm, HTi 4.0 – 5.5 mm (K072071)	IPS e.max CAD Abutment Solutions is intended for use in partially or fully edentulous mandibles and maxillae in support of single cement-retained restorations. The system comprises three parts: IPS e.max CAD mesostructure Ti base and CAD/CAM software The IPS e.max CAD mesostructure cemented to the Ti base is recommended for two-piece hybrid abutments for single tooth restorations and hybrid abutment crowns, used in conjunction with endosseous dental implants. The compatible Implant systems, Ti bases and CAD/CAM systems are shown below: Implant systems: ... Straumann Bone Level (K053088, K062129,	Equivalent



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Substantial Equivalence Discussion

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
K Number	Straumann® Variobase™ Abutments Subject Submission	IPS e.max® Press – Abutment Solutions (K120053) including P.004 RC Cementable Abutments (K072071)	IPS e.max CAD Abutment Solutions (K132209) including Sirona Dental CAD/CAM System (K111421)	
			<p>K060958) ...</p> <p>CAD/CAM Systems: Sirona inLab and Cerec SW 4.2 and above</p> <p>...</p> <p>Straumann, Bone Level NC, Ø3.3 mm, S BL 3.3, 6308154, L</p> <p>Straumann, Bone Level RC, Ø4.1 mm, S BL 4.1, 6308337, L</p> <p>...</p> <p>For the titanium base Straumann Bone Level 3.3 L the indication is restricted for replacement of single lateral incisors in the maxilla and lateral and central incisors in the mandible.</p>	
Ti-base Material	Titanium-Aluminum-Niobium alloy	Grade 4 Titanium	Titanium-Aluminum-Vanadium alloy (Ti-6Al-4V), medical grade 5, ASTM 136	Equivalent
Abutment Diameter	(b)(4) Trade Secret Process - Product Specs			Equivalent
Abutment Height	(b)(4) Trade Secret Process - Product Specs			Equivalent

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Substantial Equivalence Discussion

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
K Number	Straumann® Variobase™ Abutments Subject Submission	IPS e.max® Press – Abutment Solutions (K120053) including P.004 RC Cementable Abutments (K072071)	IPS e.max CAD Abutment Solutions (K132209) including Sirona Dental CAD/CAM System (K111421)	
Bonding Surface Area	(b)(4) Trade Secret Process - Product Specs	43 mm <sup>2</sup> (RC)	53.5 mm <sup>2</sup> (RC)	Equivalent
Coping/ Crown Material	(b)(4) Trade Secret Process - Product Specs	IPS e.max® Press Ceramic	IPS e.max® CAD Ceramic	Identical for the IPS e.max® Press Ceramic and the IPS e.max® CAD Ceramic
Design Workflow	(b)(4) Trade Secret Process - Product Specs			
Manufacturing Workflow	(b)(4) Trade Secret Process - Product Specs			
Mode of Action	(b)(4) Trade Secret Process - Product Specs			
Reusable	No	No	No	Identical

Table 10 - Comparison Matrix: Proposed Device versus Predicate Devices (K120053/K072071 and K132209/K111421)

## Traditional 510(k) Submission

### Straumann® Variobase™ Abutments

#### Substantial Equivalence Discussion

The table below provides a comparison matrix of the proposed and predicate devices (K072055):

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
<b>K Number</b>	<b>Straumann® Variobase™ Abutments Subject Submission</b>	<b>Lava™ Frame, Lava™ Frame Shade (K072055)</b>	
<b>Indications for Use</b>	The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth restorations.	<p>The Lava™ system is intended for CAD/CAM fabrication of all-ceramic dental restorations. The system is used for the manufacturing of inlays, onlays, veneers, crowns and bridges.</p> <p>Lava™ Frame and Lava™ Frame Shade are intended for the manufacturing of abutments. The titanium connection for the abutment must meet the following dimensions:</p> <ul style="list-style-type: none"> <li>- Overall cementation surface &gt;30 mm<sup>2</sup></li> <li>- Height of the head of the titanium interface from the shoulder &gt;2.8 mm</li> </ul> <p>The following systems fulfill the above described specifications:</p> <p>...</p>	Equivalent
<b>Ti-base Material</b>	Titanium-Aluminum-Niobium alloy	Titanium	Equivalent
<b>Abutment Diameter</b>	(b)(4) Trade Secret Process - Product Specs		Equivalent
<b>Abutment Height</b>	(b)(4) Trade Secret Process - Product Specs		Equivalent

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
 Substantial Equivalence Discussion

FEATURE	PROPOSED DEVICE	PREDICATE DEVICE	EQUIVALENCE DISCUSSION
K Number	Straumann® Variobase™ Abutments Subject Submission	Lava™ Frame, Lava™ Frame Shade (K072055)	
Bonding Surface Area	(b)(4) Trade Secret Process - Product Specs	≥33 mm <sup>2</sup>	Equivalent
Coping/ Crown Material	(b)(4) Trade Secret Process - Product Specs	Lava™ Frame and Lava™ Frame Shade	Identical for the Lava™ Frame and Lava™ Frame Shade
Design Workflow	(b)(4) Trade Secret Process - Product Specs		
Manufacturing Workflow	(b)(4) Trade Secret Process - Product Specs		
Mode of Action	Screw-retained or cement retained	Cement retained	Equivalent
Reusable	No	No	Identical

**Table 11 - Comparison Matrix: Proposed Device versus Predicate Devices (K072055)**

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Proposed Labeling

---

## 12 Proposed Labeling

### 12.1 Package Label

There are no changes to the package label as a result of the proposed change in this premarket notification. To aid in the review of the submission, an example of the label is shown in Figure 1.





Figure 1 - Example label for Straumann® Variobase™ Abutment

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

Proposed Labeling

---

### 12.2 Proposed Package Insert/Instructions for Use for Straumann® Variobase™ Abutments

Instructions for use: Straumann® Variobase™ Abutments




Institut Straumann AG, Peter-Merian-Weg 12, CH-4002 Basel/Switzerland,

[www.straumann.com](http://www.straumann.com)

**English**

**CAUTION: Federal law restricts this device to sale by or on the order of a dental professional.**



#### 1. Product Description

##### Abutments

Abutments are placed into dental implants to provide support for prosthetic reconstructions such as crowns.

##### Basal Screws

Basal screws are used for the fixation of the abutment to the dental implant.

#### 2. Intended use

Prosthetic components directly connected to the endosseous dental implant are intended for use as an aid in prosthetic rehabilitations.

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Proposed Labeling

---

**3. Indications**

The Straumann® Variobase™ Abutment is a titanium base placed onto Straumann dental implants to provide support for customized prosthetic restorations. Straumann® Variobase™ Abutments are indicated for screw-retained single tooth or cement-retained single tooth restorations.

**4. Contraindications**

Allergies or hypersensitivity to materials used, which are indicated in the following table:

Straumann Variobase™ Abutments, Screws	Titanium alloy, Ti-6Al-7Nb (titanium-aluminum-niobium or TAN).

**5. Warnings and Precaution**

Our products must be secured against aspiration when used intraorally. Failure to follow the procedures outlined in these instructions may lead to any or all of the following complications:

- Aspiration or swallowing of a component
- Breakage
- Infection

The Straumann® Variobase™ Abutments are single use devices.

Place implant-borne restorations only in occlusion when the implant is completely osseointegrated.

Angled abutments should not be used in areas of high mechanical loads on small diameter implants.

Dental cement or any other material used for the attachment of prosthetic components should be processed as specified by the manufacturer.

Implants are only to be restored with the corresponding original abutment compatible with that specific implant.

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Proposed Labeling

The Straumann® Variobase™ Abutments have not been evaluated for safety and compatibility in the magnetic resonance environment. The Straumann® Variobase™ Abutments have not been tested for heating or migration in the magnetic resonance environment.

#### 6. Compatibility information

Straumann® implants and the prosthetic components are available in a variety of configurations to meet your clinical needs. The label on each product uses abbreviations to help you identify whether a particular abutment or coping is compatible with the implant that you are restoring. The implant as well as the prosthetic component contains an identifier for the connection type, summarized in the table below.

Implant connection type	Compatible prostheses
NC (Narrow CrossFit®)	parts labeled NC
RC (Regular CrossFit®)	parts labeled RC
NNC (Narrow Neck CrossFit®)	parts labeled NNC
RN (Regular Neck)	parts labeled RN
WN (Wide Neck)	parts labeled WN

#### 7. Cleaning and Disinfection

Straumann® Variobase™ Abutments and components are non-sterile when delivered. Before placing the restoration in the patient's mouth, the product must be cleaned, disinfected and sterilized. Straumann recommends the following procedure for cleaning, disinfection and sterilization of abutments prior to use.

- 1) Clean rinsing under flowing water while brushing outer and inner side with adequate brushes.
- 2) The pre-treated product is to be cleaned/disinfected in an automated washer disinfector. Select the appropriate program according to the manufacturer's instructions.

#### 8. Sterilization

The restoration may be sterilized unwrapped or can be placed in an accessory cassette and packaged twice in common sterilization wraps (paper/film bags). Steam sterilize according to the parameters below:



# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Proposed Labeling

Material	Method	Conditions
Variobase™ Abutment, TAN Screw, TAN	Autoclave (moist heat) Displacement: gravity or fractionated vacuum	134 °C (273 °F) 5 minutes

**Please note:** User should ensure the use of the appropriate biological indicator for the sterilizer and parameters used.

**Please note:** User should consult the coping/restoration material manufacturer's recommendations regarding sterilization.

**Caution:** Use devices immediately after sterilization. Do not store sterilized devices.

### 9. Procedure

#### Use and handling of the Straumann® Variobase™ Abutments for the Dental Technician

##### Restoration design

When using pressing or casting techniques via wax-up, use the burn-out coping for Straumann® Variobase™ Abutments which supports a clean and sharp-edged finish of the screw channel and a good fit to the Straumann® Variobase™ Abutment. When using a digital workflow, use the Straumann Variobase Implant Kit with any CAD software platform, to facilitate the precise design of the interface between the Straumann® Variobase™ Abutment and the coping. The kit consists of an STL file containing the required milling template for the inner coping geometry. Once the coping is designed using the CAD software, send to Straumann for milling.

**Please note:** For in-lab casting, Straumann recommends using materials Type 4 dental metals (ISO 22674). The  $R_{p0.2}$  must be  $\geq 360$  MPa and have a minimum of 2% elongation after fracture. The manufacturing of the coping and/or crown must follow the standard procedure according to the material manufacturer's instructions for use.

**Please note:** For in-lab pressing, Straumann recommends using the IPS e.max® Press Abutment Solutions ceramic material according to the material manufacturer's instructions for use.

**Please note:** For milling, the following framework guidelines must be followed when designing in CAD software:

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Proposed Labeling

	Ceramic Material	Polymer Material
Minimum Wall Thickness	0.4 mm	0.5 mm
Minimum Angulation	0°	0°
Maximum Angulation	30°	30°

**Please note:** The milled restorations may be manufactured using polycon® ae, Zerion®, IPS e.max® CAD Abutment Solutions ceramic material, or Lava™ Frame. The design of the coping and/or crown must follow the standard procedure according to the material manufacturer's instructions for use.

#### Processing

Process the manufactured coping or crown following standard procedure according to the material manufacturer's instructions. Always finalize the crown or coping prior to bonding to the Straumann® Variobase™ Abutments.

#### Bonding

**Please note:** It is not necessary to sandblast the Straumann® Variobase™ Abutment.

- 1) Fix the abutment to the implant analog with a screw (hand-tight).
- 2) Seal the screw channel with wax.
- 3) Apply self-adhesive dental cement on the abutment. Only suitable self-adhesive cementation systems for the material used shall be used. Follow the instructions for use of both the dental material and cement/bonding material manufacturer. (Straumann® recommends Panavia™ F2.0 resin cement by Kuraray)
- 4) Bond the coping to the abutment.
- 5) Immediately remove excess cement from the abutment and polish the lower margin of the coping after the cement is set.
- 6) Optional: For cement retained-restorations: Make a crown following standard procedure according to the material manufacturer's instructions and finalize it.
- 7) Clean the restoration prior to sending it to the dentist.
- 8) Include this instruction for use when sending the restoration to the dentist.

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Proposed Labeling

---

#### Use and handling of the Straumann® Variobase™ Abutments for the Dentist

Remove the restoration from the master cast or the analog.

Clean, disinfect and sterilize the device as described in sections 7 and 8 of this Instructions for Use document.

#### Placing the restoration

- a) Remove the healing cap or temporary restoration.
- b) Clean and dry the interior of the implant and the abutment thoroughly.
- c) Place the sterilized restoration into the patient's mouth.
- d) Make sure that the retentive elements of the implant abutment connection are properly aligned.
- e) Use the screw delivered with the abutment to screw the abutment into the dental implant.

**Please note:** Always ensure that the surfaces of threads and screw heads are clean and that a new screw is used for the restoration.

- f) Straumann® abutments are fixed to the implant using the Straumann® SCS screwdriver, ratchet and torque control device. Use the respective torque according to the table below:

Device type	Tightening torque	Special considerations
Abutments (permanent)	35 Ncm	n/a
Temporary abutments	15 – 35 Ncm	Tighten only to 35 Ncm if the implant is fully osseointegrated
Components on implant analogs	Hand-tight	n/a

For cement-retained restorations (optional):

- g) Close the screw channel with cotton and sealing compound (i.e., gutta-percha)
- h) Apply self-adhesive dental cement on the two-piece abutment. Only suitable self-adhesive cementation systems for the used materials shall be used. Follow the instructions for use of the cement/bonding material manufacturer (Straumann recommends Panavia™ F2.0 resin cement by Kuraray).

# Traditional 510(k) Submission

## Straumann® Variobase™ Abutments

### Proposed Labeling

---

- i) Bond the crown to the two-piece abutment.
- j) Immediately remove excess cement from the two-piece abutment.

#### **Warning**

Torques greater than 35 Ncm may result in failure of the abutment and/or implant. Torque values less than the recommended values may result in loosening of the abutment, which may lead to abutment and/or implant failure.

#### **Please note**

Once the Straumann® abutment has been secured to the implant using the indicated torque, it should not be removed.

#### **10. Further Information**

For additional information about the use of Straumann® products, call Straumann's customer service department or visit [www.straumann.com](http://www.straumann.com).

For additional information, consult:

Basic information on the Straumann® Variobase™ Abutment

#### **11. Please note**

Practitioners must have appropriate knowledge and instruction in the handling of the Straumann product described herein ("Straumann Product") for using the Straumann Product safely and properly in accordance with these instructions for use.

The Straumann Product must be used in accordance with the instructions for use provided by the manufacturer. It is the practitioner's responsibility to use the device in accordance with these instructions for use and to determine, if the device fits to the individual patient situation.

The Straumann Product is part of an overall concept and must be used only in conjunction with the corresponding original components and instruments distributed by Institut Straumann AG, its ultimate parent company and all affiliates or subsidiaries of such parent company ("Straumann"), except if stated otherwise in these instructions for use. If use of products made by third parties is not recommended by Straumann in these instructions for use, any such use will void any warranty or other obligation, express or implied, of Straumann.

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**  
Proposed Labeling

---

**12. Validity**

Upon publication of these instructions for use, all previous versions are superseded.

© Institut Straumann AG, 2013. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

Panavia™ is a trademark of Kurary Co, LTD, JP.

**Availability**

Some items of the Straumann® Dental Implant System are not available in all countries.



Straumann Products with the CE mark fulfill the requirements of the Medical Devices Directive 93/42 EEC /



Refer to package insert



Manufacturer



Article number



Lot Number



Do not re-use



Non-sterile

Rx only

Federal law restricts this device to sale by or on the order of a dental professional.


**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 2

---

**Appendix 2 – Basic Information on the Straumann® Variobase™  
Abutment**









(b)(4) Trade Secret Process - Draft Product Specs





(b)(4) Trade Secret Process - Draft Product Specs





(b)(4) Trade Secret Process - Draft Product Specs





(b)(4) Trade Secret Process - Draft Product Specs





(b)(4) Trade Secret Process - Draft Product Specs



(b)(4) Trade Secret Process - Draft Product Specs





(b)(4) Trade Secret Process - Draft Product Specs





(b)(4) Trade Secret Process - Draft Product Specs





(b)(4) Trade Secret Process - Draft Product Specs







(b)(4) Trade Secret Process - Draft Product Specs



**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**


Appendix 16

---

**Appendix 16 – IPS e.max® Press Abutment Solutions**  
**Instructions for Use**

# X<sup>®</sup>Press

Abutment Solutions



# e IPS



all ceramic  
all you need

## INSTRUCTIONS FOR USE

CE 0123

Traditional 510(k)  
Straumann® Variobase™ Abutments

Straumann USA, LLC  
January 23, 2014

ivoclar  
vivadent  
technical

A16-2

# Table of Contents

	3	<b>Product Information</b> Material Uses Composition
APPLICATION PROCEDURE	6	<b>Fabrication of Hybrid Abutments and Hybrid Abutment Crowns</b> Treatment/fabrication process Shade – tooth shade, preparation shade/abutment shade Model preparation Selecting a Ti base Layer thicknesses of the ceramic components Modeling Sprueing Investing Preheating Pressing Divesting Removing the reaction layer Finishing Stain and Characterization firing Glaze firing Crown on hybrid abutment
	32	<b>Optional: Clinical Try-in</b> Temporarily securing the pressed object on the Ti base Clinical try-in
	37	<b>Permanent Cementation</b> Pre-treatment of the Ti base Preparing the pressed object Cementation with Multilink Implant
	43	<b>Seating and Aftercare</b> Sterilization Intraoral preparation Seating the hybrid abutment and crown Seating the hybrid abutment crown Care notes – Implant Care
	50	<b>General Information</b> Frequently Asked Questions Material Selection Table Press and Firing Parameters Clinical Cases

# IPS e.max<sup>®</sup> Press Abutment Solutions

## Product Information

Press ceramics have been synonymous with esthetics, accuracy of fit, shape and function for decades. The IPS e.max Press lithium disilicate (LS<sub>2</sub>) glass-ceramic additionally offers an outstanding strength of 400 MPa. The already extensive indication range from thin veneers (0.3 mm) and monolithic molar crowns to anterior and premolar bridges is now expanded to include hybrid abutment restorations.

With IPS e.max Press, you can fabricate such restorations in combination with a titanium base (Ti base). Two different approaches are available:

- Hybrid abutment and separate crown
- Hybrid abutment crown

Both solutions show outstanding function, efficiency and esthetics. The durable bond to the Ti base is achieved by means of the self-curing Multilink Implant luting composite.

### Hybrid abutment

The hybrid abutment is an individually pressed LS<sub>2</sub> abutment which is luted to the Ti base. The shape, emergence profile and esthetic properties of this abutment can be ideally adjusted to the clinical situation.

### Individual esthetics

Given the lifelike appearance of LS<sub>2</sub> glass-ceramics, the esthetic possibilities are virtually limitless, particularly in the anterior region. Due to the individual characterization, a lifelike appearance is achieved near the root and the transition area to the crown. With the preparation margin of the crown located on the gingival level, the geometry of the hybrid abutments allows for an easy integration of the restoration. Excess cementation material is therefore easily removed.

### Flexibility due to laboratory fabrication

The pressed LS<sub>2</sub> abutment is extraorally luted to a Ti base with Multilink Implant, then screwed into place in the oral cavity and finally provided with a permanent IPS e.max crown. As the hybrid abutment is conveniently fabricated in the lab, the process is time-saving and flexible.

### Hybrid abutment crown

Hybrid abutment crowns are characterized by combining abutment and monolithic crown in one piece. This is an efficient two-in-one solution made of lithium disilicate (LS<sub>2</sub>), which is directly luted to a Ti base.

### Efficient fabrication due to two-in-one approach

LS<sub>2</sub> glass-ceramics provide for strength, durability and efficiency, particularly in the posterior area. Moreover, the material offers well-known esthetic properties allowing restorations to be simply characterized with IPS e.max Ceram stains.

### Luted extraorally, screwed in intraorally – for improved flexibility


The monolithically pressed hybrid abutment crown is reliably luted to the Ti base by means of Multilink Implant. Then, the restoration is screwed onto the implant – in one piece. Thus, the bothersome task of excess cement removal is a thing of the past. Subsequently, the screw access channel is sealed with a composite material (e.g. Tetric EvoCeram<sup>®</sup>). If required, the screw can be accessed at any time, which affords the dental team clinical flexibility.

### New possibilities for economically efficient restorations

IPS e.max Press hybrid abutment crowns are a new, economically attractive alternative to conventional implant-supported restorations, particularly for the posterior area, where strength, durability and convenient clinical handling matter.

### Note regarding the Instructions for Use:

The present Instructions for Use deal only with IPS e.max Press Abutment Solutions and represent a supplement to the existing IPS e.max Press Instructions for Use. The IPS e.max Press Instructions for Use contain more detailed descriptions of the material (e.g. ingot concept) and the entire indication range. In case you do not have the IPS e.max Press Instructions for Use, you can order them from your sales representative or simply download it from [www.ivoclarvivadent.com](http://www.ivoclarvivadent.com).




## Material

### Press ceramic

IPS e.max Press are lithium disilicate glass-ceramic (LS<sub>2</sub>) ingots for the press technology. The industrial production process creates absolutely homogeneous ingots in different translucency levels. For IPS e.max Press Abutment Solutions, ingots from the existing range of products are used.

The ingots feature a strength of 400 MPa and are thus the pressed ceramic ingots with the highest strength. They are pressed in Ivoclar Vivadent press furnaces to produce objects with outstanding accuracy of fit. The pressed, tooth-coloured, highly esthetic restorations are completed with IPS e.max Ceram.



**IPS e.max Press** Lithium-Disilicate

CTE (100–400°C) [10 <sup>-6</sup> /K]	10.2
CTE (100–500°C) [10 <sup>-6</sup> /K]	10.5
Flexural strength (biaxial) [MPa]*	400
Fracture toughness [MPa m <sup>0.5</sup> ]	2.75
Modulus of elasticity [GPa]	95
Vickers hardness [MPa]	5800
Chem. solubility [μg/cm <sup>2</sup> ]*	40
Pressing temperature [°C]	915–920

\*according to ISO 6872

### Ti base

For IPS e.max Press Abutment Solutions, customary Ti bases made of titanium or titanium alloys are used.

Please observe the instructions for use and processing of the manufacturer of the Ti bases used.

## Uses

### Indications

- Hybrid abutments for anterior and posterior single-tooth restorations
- Hybrid abutment crowns for anterior and posterior restorations

### Contraindications

- Use of Ti bases which do not fulfil the geometry requirements.
- Failure to observe the requirements stipulated by the implant manufacturer for using the selected implant type (diameter and length of the implant must be approved for the respective position in the jaw by the implant manufacturer).
- Bruxism
- Failure to observe the permissible maximum and minimum ceramic wall thicknesses.
- All uses not stated as indications are contraindicated.

### Important processing restrictions

Failure to observe the following restrictions may compromise the results achieved with IPS e.max Press:

- If hybrid abutment crowns are fabricated, the opening of the screw channel must not be located in the area of contact points and areas with masticatory function. If this is not possible, a hybrid abutment with a separate crown would be preferred.
- No extension units; only single-tooth restorations
- Layering with a veneering ceramic other than IPS e.max Ceram
- Pressing of two or more IPS e.max Press ingots in one investment ring
- Pressing of IPS e.max Press in the IPS Investment System 300 g
- Use of a luting composite other than Multilink® Implant to lute IPS e.max Press to the Ti base
- Temporary cementation of the crown on the hybrid abutment.
- Failure to observe the manufacturer's instructions regarding the processing of the Ti base.

### Side effects

If the patient is known to be allergic to any of the components, IPS e.max Press Abutment Solutions should not be used.

## Composition

IPS e.max Press ingots and accessories required in conjunction with IPS e.max Press Abutment Solutions consist of the following main components:

- **IPS e.max Press ingots**  
Components: SiO<sub>2</sub>  
Additional components: Li<sub>2</sub>O, K<sub>2</sub>O, MgO, ZnO, Al<sub>2</sub>O<sub>3</sub>, P<sub>2</sub>O<sub>5</sub> and other oxides
- **IPS AloX Plunger**  
Components: Al<sub>2</sub>O<sub>3</sub>
- **IPS AloX Plunger Separator**  
Components: Boron nitride
- **IPS e.max Press Invex Liquid**  
Components: Hydrofluoric acid and sulphuric acid in water
- **IPS PressVEST Powder**  
Components: SiO<sub>2</sub>, MgO and NH<sub>4</sub>H<sub>2</sub>PO<sub>4</sub>
- **IPS PressVEST Liquid**  
Components: Colloidal silicic acid in water
- **IPS PressVEST Speed Powder**  
Components: SiO<sub>2</sub>, MgO and NH<sub>4</sub>H<sub>2</sub>PO<sub>4</sub>
- **IPS PressVEST Speed Liquid**  
Components: Colloidal silicic acid in water
- **IPS Object Fix Flow**  
Components: Oxides, water, thickening agent
- **IPS Ceramic Etching Gel**  
Components: Hydrofluoric acid (approx. 5%)
- **Virtual Extra Light Body Fast Set**  
Components: Addition-reaction silicone, vinyl polysiloxane, methylhydrogensiloxane, organoplatinic complex, silica
- **Monobond Plus**  
Components: Alcohol solution of silane methacrylate, phosphoric acid methacrylate and sulphide methacrylate
- **Multilink Implant**  
Components: Dimethacrylate, HEMA, barium glass, ytterbium tri-fluoride, spheroid mixed oxide


## Warning

- Do not inhale ceramic dust during finishing – use exhaust air discharge and mouth protection.
- IPS Ceramic Etching Gel contains hydrofluoric acid. Contact with skin, eyes and clothing must be prevented at all costs, since the material is extremely toxic and corrosive. The etching gel is intended for professional use only and must not be applied intraorally (inside the mouth).



# IPS e.max<sup>®</sup> Press Abutment Solutions

## Treatment/Fabrication Process



## Shade – tooth shade, preparation shade/abutment shade

Optimum integration in the oral cavity of the patient is the prerequisite for a true-to-nature all-ceramic restoration. To achieve this, the following guidelines and notes must be observed by both the dentist and the laboratory.

With IPS e.max Press Abutment Solutions, you can imitate not only the clinical crown of a natural tooth, but also a part of the "root". By defining/determining the "root shade" you can adjust the shade of the IPS e.max Press Abutment Solution restoration accordingly. This allows you to achieve a highly esthetic implant-supported restoration which retains its lifelike appearance also in the case of gingiva recession.

Hybrid abutment and separate crown		
<b>Restoration shade</b> (pressed ceramic LS <sub>2</sub> , characterization)	<b>Shade luting material</b> (crown on hybrid abutment)	<b>Shade hybrid abutment</b> (Ti base, luting material, pressed ceramic LS <sub>2</sub> )
		
Hybrid abutment crown		
<b>Shade hybrid abutment crown</b> (Ti base, luting material, pressed ceramic LS <sub>2</sub> , characterizations)		

Please refer to the table on page 52 for the selection of the suitable IPS e.max Press ingot.

## Model preparation

For the fabrication of an IPS e.max Press Abutment Solutions restoration, a model with gingiva mask is fabricated.

- Select the suitable model analog according to the implant system used.
- Fabricate a model with gingiva mask.




Prepared model with gingiva mask

## Selecting a Ti base

The following paragraphs outline the selection criteria for a suitable Ti base. As a general rule, the instructions of the respective manufacturer regarding the use of the Ti base have to be observed.

- Only bases consisting of Ti or Ti alloys must be used.
- Select a Ti base with a size that matches the clinical situation and the chosen implant system. The geometry requirements must be observed.
- The rotation lock must be designed in such a way that stress concentrations on the pressed object are avoided.
- Ti bases with undercuts, e.g. retention grooves, are suitable to some extent.
- Check the available space for the pressed object taking the geometry of the Ti base into account on the model (e.g. silicone key).
- Observe the instructions of the manufacturer when modifying the Ti base.

	Minimum dimensions	
	Height $H_{Ti}$ (bonding surface)	Shoulder width $S_{Ti}$
Ti base	$H_{Ti}$ min. 4.0 mm	$S_{Ti}$ min. 0.6 mm



## Layer thicknesses of the ceramic components

Observing the geometry requirements of the pressed objects made of IPS e.max Press material is the key to success for a durable restoration. The more attention given to the design, the better the final results and the clinical success will turn out to be. The following basic guidelines have to be observed:

### Hybrid abutment

- The marginal shoulder width  $S_A$  must be at least 0.6 mm.
- Create an emergence profile with a right angle at the transition to the crown (see picture).
- The wall thickness  $W_A$  must be at least 0.5 mm.
- The height  $H_A$  must not exceed twice the height of the Ti base  $H_{Ti}$ .
- The hybrid abutment should be designed in a similar way as a prepared natural tooth:
  - Circular epi-/supragingival shoulder with rounded inner edges or a chamfer.
  - In order for the crown to be cemented to the hybrid abutment using a conventional/self-adhesive cementation protocol, retentive surfaces and a sufficient "preparation height" must be observed.
- The width  $B_{AK}$  of the crown is limited to 6.0 mm from the axial height of contour to the screw channel of the hybrid abutment.

### Hybrid abutment crown

- The marginal shoulder width  $S_A$  must be at least 0.6 mm.
- The wall thickness  $W_{AK}$  must be larger than 1.5 mm for the entire circumference.
- The opening of the screw channel must not be located in the contact point areas or areas with a masticatory function. If this is not possible, a hybrid abutment with a separate crown should be preferred.
- The width of the hybrid abutment crown  $B_{AK}$  is limited to 6.0 mm from the axial height of contour to the screw channel.
- The height  $H_{AK}$  must not exceed twice the height of the Ti base by more than 2 mm.

## Modeling

### Fabrication of a resin coping

To prepare the wax-up, a resin coping is prepared if both hybrid abutments and hybrid abutment crowns are fabricated. Please observe the following procedure:

- Check the implant position and inclination with regard to the position of the screw channel.
- Screw the Ti base onto the model analog with the corresponding screw.
- **Tip:** Make sure that an additional model analog is available as this will facilitate some steps.
- Clean the Ti base with a steam cleaner.
- Insert a pin with the same diameter as the screw channel to "seal" and "extend" the screw channel.
- Do not apply die spacer.
- Isolate the Ti base and the pin with a thin application of separator. If too much separator is used, this might result in uneven areas on the inner aspect of the pressed object.
- In order to achieve a sound fit and to facilitate the subsequent wax-up, a coping is first fabricated on the Ti base with modelling resin. Design the coping in such a way that it can subsequently be completely covered with modelling wax. Please observe the instructions of the manufacturer regarding the processing of modelling resin.
- Remove the Ti base from the model.
- Eliminate possible over-contoured areas of the resin coping at the transition area to the Ti base by means of rubber polishers. Do not damage the Ti base.
- Remove the resin coping together with the pin from the Ti base.
- Loosen and remove the pin by rotating the resin coping.
- Screw the Ti base onto the model analog again.
- Place the resin coping back on the Ti base and check the fit and dimension (e.g. silicone key). If necessary, adjust the coping by means of rotary instruments.



Screw the Ti base onto the model analog with the corresponding screw.



- Insert a pin with the same diameter as the screw channel to "seal" and "extend" the screw channel.



Isolate the Ti base and the pin with a thin application of separator.




Apply the modelling resin to the Ti base in increments.




Design the resin coping on the entire Ti base.



Remove the resin coping together with the pin from the Ti base.




Loosen and remove the pin by rotating the resin coping.



Eliminate possible over-contoured areas of the resin coping at the transition area to the Ti base by means of rubber polishers.



Place the resin coping back on the Ti base and check the fit and dimension (e.g. silicone key). If necessary, adjust the coping by means of rotary instruments. Design the coping in such a way that it can subsequently be covered with modelling wax.



### Wax-up

Please observe the following notes with regard to modelling:

- Observe the stipulated layer thicknesses.
- Create an accurate model of the restoration, particularly at the transition area to the Ti base.
- Do not over-contour the margins, since this would require time-consuming and risky fitting procedures after pressing.
- Use an organic wax for modelling to ensure that it burns out without leaving residue in the investment ring.

### Procedure for hybrid abutments

- Before creating the wax object, re-insert the isolated pin into the screw channel.
- Design the emergence profile by flooding the area between the gingiva mask and the resin coping with wax.
- Contour the hybrid abutment to a reduced tooth shape. The hybrid abutment should be designed in such a way that the required layer thicknesses are met in the crown that is fabricated. Check by means of the silicone key and in relation to the opposing dentition.
- Determine the crown margins in relation to the gingiva level.
- Design a chamfer on which the crown is subsequently seated.
- Remove the object together with the Ti base from the model and check the emergence profile. If necessary, make adjustments.
- Check the transition to the Ti base and remove excess wax.
- Check the required minimum thicknesses (page 9) prior to attaching the sprue.



Design the emergence profile by flooding the area between the gingiva mask and the resin coping with wax.



Design the hybrid abutment with a reduced tooth shape and determine the crown margin in relation to the gingiva level.



Check the dimensions by means of the silicone key and in relation to the opposing dentition.



Remove the object together with the Ti base from the model and check the emergence profile. If necessary, make adjustments. Check the transition to the Ti base and thoroughly remove excess wax.



**Procedure for hybrid abutment crowns:**

- If required, re-insert the isolated pin into the screw channel before creating the wax object.
- Design the emergence profile by flooding the area between the gingiva mask and the resin coping with wax.
- Design the abutment crown to full contour according to functional and esthetic criteria. Check in relation to the opposing dentition.
- Make sure to take a slightly reduced occlusal relief into consideration during the wax-up, since the application of the Stains and Glaze results in a slight increase in vertical dimensions.
- Remove the object together with the Ti base from the model and check the emergence profile. If necessary, make adjustments.
- Check the transition to the Ti base and remove possible excess wax.
- Check the required thicknesses (page 9) prior to attaching the sprues.



Design the emergence profile by flooding the area between the gingiva mask and the resin coping with wax. If necessary, re-insert the pin prior to modelling.




Design the abutment crown to full contour according to functional and esthetic criteria. Check the object in relation to the opposing dentition.




## Sprueing

Please observe the following notes when attaching the sprues to the abutment or the abutment crown:



- Depending on the number and size of the objects to be invested, either the 100 g or 200 g IPS Investment Ring System is selected. Before sprueing, weigh the ring base and record the weight (seal the opening of the ring base with wax). Please note that the mixing ratio of the investment material is different for the various restoration types (e.g. inlays, crowns, abutments).
- Use a **2.5 mm wax wire** for sprueing.
- For abutments, the sprue is attached to an axial surface.
- For abutment crowns, the sprue is attached to a cusp.
- **Align the wax wire as parallel as possible to the screw channel in order to prevent the investment material from fracturing in the screw channel.**
- The maximum length (object + sprue) of 16 mm must not be exceeded.
- Place the object on the investment ring base in such a way that the screw channel is parallel to the outer wall of the investment ring. As a result, the investment material can subsequently be filled evenly and in a controlled manner. The objects could be placed in a tilted position on the investment ring base, but this may lead to difficulties during investing (e.g. bubbles in the screw channel).
- Observe a distance of at least 10 mm between the object and the silicone ring.
- If only one object is invested and pressed in an EP500 furnace, a second short (blind) sprue must be placed. This ensures that the switch-off function of the furnace works properly at the end of the pressing procedure.



Attach the sprue to a circular area of the abutment model and as parallel as possible to the screw channel. Use a 2.5 mm wax wire.



Attach the sprue to an oral cusp of the abutment crown and as parallel as possible to the screw channel. Use a 2.5 mm wax wire.



Place the object on the investment ring base in such a way that the screw channel is parallel to the outer wall of the investment ring. As a result, the investment material can subsequently be filled evenly and in a controlled manner. Furthermore, this reduces the risk of the investment material breaking in the screw channel while the ceramic material is pressed.

## Investing

Investing is carried out with either IPS PressVEST or IPS PressVEST Speed. The corresponding IPS Silicone Ring with the matching ring gauge is used for investment.

Determine the weight of the object before investing.

- Position the wax objects on the ring base and attach them with wax and weigh.
- The difference between the empty and the loaded ring base is the definitive wax weight.

	Small Ingot	Large Ingot (L)
Wax weight	up to max. 0.75 g	up to max. 2 g
Investment Ring System	100 g and 200 g	only 200 g

Please refer to the Instructions for Use of the corresponding investment material regarding the detailed processing parameters. The following basic procedure is recommended:

- Do not use a debubbler on the wax objects.
- Carefully place the IPS Silicone Ring on the ring base without damaging the objects. The silicone ring must sit flush on the ring base.
- The processing temperature of the investment material is 18– max. 23 °C / 64 °F – max. 73 °F. Higher or lower processing temperatures substantially affect the setting behaviour.
- Mix the investment material. Note: The investment material contains quartz powder. Therefore, avoid the inhalation of dust.
- **Important: Pour the investment material slowly into the investment ring, so that the material continuously fills the screw channel. If the material does not sufficiently fill the screw channel, use an instrument to carefully apply additional investment material to the screw channel from the top.**
- Carefully fill the investment ring with investment material up to the marking on the silicone ring and position the ring gauge with a hinged movement.
- Allow the investment ring to set without manipulating it.
- To prevent crystallization of the IPS PressVEST investment material, the invested ring must be processed within 24 hours.
- If IPS PressVEST Speed is used, make sure that the investment ring is placed in the preheating furnace after a setting time of at least 30 and maximum 45 minutes.

### Investment material: Liquid concentration and quantity

Indication	IPS PressVEST		IPS PressVEST Speed	
	100 g Investment Ring Liquid : dist. water	200 g Investment Ring Liquid : dist. water	100 g Investment Ring Liquid : dist. water	200 g Investment Ring Liquid : dist. water
<b>IPS e.max Press</b>				
Hybrid abutment Hybrid abutment crown	16 ml : 6 ml	32 ml : 12 ml	20 ml : 7 ml	40 ml : 14 ml
<b>Mixing time</b> (under vacuum at approx. 350 rpm)	60 seconds		2.5 minutes If a high-speed mixer is used, the mixing time under vacuum has to be reduced.	

**Liquid concentration:** The data contained in the table are approximative values. Depending on the geometry of the Ti base and the materials used for the wax-up, these values may be individually changed. However, the concentrated Liquid content must not be lower than 50% in relation to distilled water.

**Important:** The total quantity of liquid (Liquid + dist. water) must not be altered.



Correctly sprued abutment (left) and abutment crown (right). The screw channel is in a vertical position and parallel with the wall of the investment ring.



Pour the investment material slowly into the investment ring, so that the material can continuously fill the screw channel.



Continue to carefully fill the investment ring up to the marking and position the ring gauge with a hinged movement.



## Preheating

After the stipulated setting time of the respective investment material (IPS PressVEST or IPS PressVEST Speed), the investment ring is prepared for preheating as follows:

- Remove the ring gauge with a turning movement.
- Carefully push the investment ring out of the IPS Silicone Ring.
- Remove the ring base with a turning movement.
- Remove rough spots on the bottom surface of the investment ring with a plaster knife. Check the 90° angle. Investment material residue must not enter the sprues. Blow into the sprues if necessary.
- If several investment rings are preheated together, mark them accordingly.

	IPS PressVEST	IPS PressVEST Speed
Setting time	min. 60 min, max. 24 hrs	min. 30 min, max. 45 min
Temperature of the preheating furnace when placing the investment ring	Room temperature	850 °C / 1562 °F; switch on the preheating furnace in time.
Position of the investment ring in the preheating furnace	Towards the rear wall, tipped with the opening facing down	Towards the rear wall, tipped with the opening facing down
Final temperature for preheating the investment ring	850 °C / 1562 °F	850 °C / 1562 °F
Holding time of the investment ring at final temperature	min. 60 min	100 g investment ring – min. 45 min 200 g investment ring – min. 60 min
IPS e.max Press ingots	<b>no preheating</b>	<b>no preheating</b>
IPS Alox Plunger	<b>no preheating</b>	<b>no preheating</b>
Important		If several Speed investments are to be conducted, they should be invested consecutively and placed into the preheating furnace at an interval of approx. 20 minutes. Make sure that the furnace temperature does not drop too much when placing the investment rings into the preheating furnace. The stipulated holding time counts from the point when the preheating temperature has been reached again.



Towards the rear wall, tipped with the opening facing down



Do not preheat the IPS e.max Press ingot and Alox Plunger.

In order to ensure smooth working procedures in the laboratory on a daily basis, impeccable functioning of the infrastructure, particularly the preheating furnaces, is essential. This includes their maintenance, cleaning with a vacuum cleaner in a cool state as well as regular checks of the temperature controls and heating elements, etc., by the manufacturer.

## Pressing

Carry out the following preparatory steps for pressing before the preheating cycle for the investment ring has been completed:

- Provide a **cold** IPS Alox Plunger and a **cold** IPS e.max Press ingot in the desired shade (please refer to the material selection table on page 52).
- Dip the **cold** IPS Alox Plunger into the opening of the IPS Alox Plunger Separator and keep it ready for use.
- Turn on the press furnace (e.g. Programat EP 5000) in time so that the self-test and preheating phase are completed.
- Select the press program for IPS e.max Press and the desired investment ring size.

Remove the investment ring from the preheating furnace immediately after completion of the preheating cycle. This step may take max. 30 seconds to prevent the investment ring from cooling down too much.

- Place the **cold** IPS e.max Press ingot into the **hot** investment ring.
- Insert the ingots in the investment ring with the non-imprinted side facing down. The imprinted side faces up to check the ingot shade.
- Place the side of the **cold** IPS Alox Plunger which has been coated with Separator into the **hot** investment ring.
- Use the investment ring tongs to place the loaded investment ring in the centre of the **hot** press furnace.
- The selected press program is started by pressing START.

After the end of the press cycle (optical and/or acoustic signal) proceed as follows:

- Remove the investment ring from the press furnace using the investment ring tongs immediately after pressing.
- Place the investment ring on a cooling grid to cool in a place protected from draft.
- Do not speed up cooling, e.g. by blasting with compressed air.

	100 g Investment Ring	200 g Investment Ring
	1 small ingot	1 small ingot or 1 large ingot
<b>IPS e.max Press ingots</b>	<b>cold ingot</b>	<b>cold ingot</b>
<b>IPS Alox Plunger</b>	<b>cold plunger</b>	<b>cold plunger</b>
<b>IPS Alox Plunger Separator</b>	✓	✓


**Ingots: select one large or one small ingot according to the determined wax weight!**




Provide a **cold** isolated IPS Alox Plunger and a **cold** IPS e.max Press ingot in the desired shade.



Place the **cold** IPS e.max Press ingot into the **hot** investment ring with the imprinted side facing up.



Then, place the side of the IPS AloX Plunger which has been coated with Separator into the **hot** investment ring.



Place the **hot** and loaded investment ring in the centre of the **hot** press furnace using the IPS Investment Ring Tongs.



Press **START** to start the selected program.



Once the press program is completed, place the hot investment ring on the cooling grid using the Investment Ring Tongs and allow it to cool to room temperature.

**Press parameters for IPS e.max Press**

Programat EP 3000



Select the press program in accordance with the selected ingot to be pressed and the investment ring size used.



The press parameters for HO, MO, LT and HT are integrated starting with software V 6.1.

Programat EP 5000



Select the press program in accordance with the selected ingot to be pressed and the investment ring size used.



The press parameters for HO, MO, LT and HT are integrated starting with software V 6.1.

The press parameters for older-generation Ivoclar Vivadent press furnaces are listed on page 53 under Press Parameters.



## Divesting

After cooling to room temperature (approximately 60 minutes), the investment ring may show cracks, which developed during the cooling phase (immediately around the Alox plunger). This is the result of the different CTEs of the various materials (Alox Plunger, investment material and press ingot) and does not compromise the press results.

Divest the investment ring as follows:

- Mark the length of the Alox Plunger on the cooled investment ring.
- Separate the investment ring using a separating disc. This predetermined breaking point enables reliable separation of the Alox Plunger and the ceramic material.
- Break the investment ring at the predetermined breaking point using a plaster knife.
- Always use polishing beads to divest the pressed objects (rough and fine divestment). Do not use  $Al_2O_3$ .
- Rough divestment is carried out with polishing beads at 4 bar (58 psi) pressure.
- Fine divestment is carried out with polishing beads at 2 bar (29 psi) pressure.
- Observe the blasting direction and distance to prevent damage to the object margins during divestment.
- Same as the outer surfaces, thoroughly blast the screw channel with polishing beads at 2 bar (29 psi) pressure.
- Remove possible ceramic residue from the Alox Plunger with type 100  $Al_2O_3$ .



Mark the length of the Alox Plunger.



Separate the investment ring using a separating disc and break it at the predetermined breaking point.

### Tip

Pull out the plunger with pliers from the separated segment using a rotating movement. This also removes any possible ceramic residue from the Alox plunger.





Rough divesting with polishing beads at 4 bar (58 psi) pressure until the object becomes visible.



Fine divestment of the abutment is carried out with polishing beads at 2 bar (29 psi) pressure.



Fine divestment of the abutment crown is carried out with polishing beads at 2 bar (29 psi) pressure.



Completely divested IPS e.max Press objects.



## Removing the reaction layer

After fine divestment, the reaction layer formed during the press procedure is removed using IPS e.max Press Invex Liquid. The procedure is carried out as follows:

- Pour the Invex Liquid into a plastic cup.
- Immerse the pressed object in the Invex Liquid and clean in an ultrasonic cleaner for at least 10 min and at most 30 min. Make sure that the objects are completely covered with Invex Liquid.
- Use the sieve insert to remove the restoration from the Invex Liquid and clean the object under running water and blow dry.
- Carefully remove the white reaction layer with type 100  $\text{Al}_2\text{O}_3$  at max. 1–2 bar (15–29 psi) pressure.
- Make sure that the reaction layer is completely removed, both from the screw channel and the outer side of the object (repeat the procedure, if necessary).
- If the reaction layer is not completely removed, problems may occur in the further course of the fabrication process.
- Replace the IPS e.max Invex Liquid after 20 applications or after sedimentation of the Liquid.



To remove the reaction layer, immerse the pressed objects in IPS e.max Press Invex ...



... and place in an ultrasonic cleaner for at least 10 and at most 30 minutes.



Using  $\text{Al}_2\text{O}_3$  and a pressure of at most 1–2 bar (15–29 psi) pressure, carefully remove the reaction layer from the outer side ...



... and the screw channel.

### Warning

- The Invex Liquid contains < 1% hydrofluoric acid.
- It is harmful when inhaled, swallowed and when it comes into contact with the skin. Furthermore, it is corrosive.
- Keep the container tightly sealed and store it in a well-ventilated place (acid cabinet).
- If the material comes into contact with the eyes, immediately rinse with copious amounts of water and see an ophthalmologist immediately.
- In case of accidental contact with skin, immediately wash with plenty of water.
- Use suitable protective clothing, gloves and goggles when working.
- In case of an accident or physical discomfort, see a physician immediately, and take the Invex label with you, if possible.



### Disposal

- Neutralize the Invex Liquid before disposal!
- Use the IPS Ceramic Neutralization Powder to neutralize the Invex Liquid.
- For 50 ml Invex Liquid, approx. 3–4 g of IPS Ceramic Neutralization Powder are required.
- Note: strong foam development during neutralization.
- Carefully add the neutralization powder to the Invex Liquid in small portions until foam is no longer formed; then allow a reaction time of 5 minutes.
- If larger quantities are disposed of, check the liquid with litmus paper (must show an alkaline reaction).
- After the reaction time, pour the neutralized solution into the sink, flushing it with running water.



## Finishing

It is of critical importance to use the correct grinding instruments for finishing and adjusting high-strength glass-ceramics (please refer to the Ivoclar Vivadent Flow Chart "Recommended grinding instruments for glass-ceramics"). If unsuitable grinding instruments are used, chipping of the edges and local overheating may occur.

Observe the following procedure for finishing IPS e.max Press restorations:

- Adjustment by grinding of pressed IPS e.max Press restorations should be kept to a minimum.
- Overheating of the ceramic must be avoided. Low speed and light pressure is recommended.
- Make sure that the minimum thicknesses are maintained even after the minor adjustments.

### Fitting to the Ti base

The fit of the abutment or abutment crown is checked on the Ti base before the sprue is separated.


- Before the object is fitted, the inner aspect of the object (screw channel) is checked for bubbles in the ceramic. If required, the bubbles are removed with suitable instruments.
- Carefully position the abutment or abutment crown on the Ti base. **Note:** Apply only light pressure to secure the pressed object on the Ti base in order to prevent chipping of the ceramic. Observe the position of the rotation lock.
- Possible rough spots interfering with the fit of the pressed object on the Ti base cause greyish-black markings on the screw channel. Carefully remove such markings with suitable grinding instruments. The diameter of the grinding instrument must be smaller than that of the screw channel. As an alternative to marking the rough spots, an occlusion spray can also be used.
- Carefully remove possible rough spots until an optimum fit between the Ti base and the pressed object is achieved. Repeat the procedure, if required.



Check the screw channel for bubbles ...



Carefully position the abutment or abutment crown on the Ti base.



Possible rough spots interfering with the fit on the Ti Base cause stains on the screw channel of the pressed object, ...



... which can be carefully removed by means suitable grinding tools.



After possible rough spots have been removed, an optimum fit between the hybrid abutment ...




... or the abutment crown and the Ti base is achieved.


### Finishing

Once an optimum fit between the abutment or the abutment crown and the Ti base has been achieved, please proceed as follows for the finishing steps:

- Separate the sprue using a separating disc. Avoid overheating.
- Smooth out the attachment point of the sprue. Make sure that the minimum thicknesses are maintained.
- Check the emergence profile and the fit on the model.
- In the case of abutment crowns, additionally check the occlusion and articulation. Adjust by grinding, if necessary, and create surface textures.
- To clean the abutment crown, briefly blast the outer side with  $Al_2O_3$  at 1 bar (15 psi) pressure and clean with the steam cleaner. Some blasting devices may require different pressure settings to accomplish this procedure.



Separate the sprues using a separating disc. Avoid overheating.



Smooth out the attachment point of the sprue.



Check the emergence profile and the fit on the model.



## Stain and Characterization firing

The following paragraphs will explain the steps of optional staining and characterizing with IPS e.max Ceram Shades and Essences. On abutments, only the emergence profile is characterized for the individual patient. This characterization may also take place at a later stage, i.e. when the crown is characterized.



If abutments are fabricated, only the area of the emergence profile is characterized with IPS e.max Ceram Shades and Essences.



If abutment crowns are fabricated, the entire outer surface may be individually characterized.

The Stains and Characterization firing is conducted with IPS e.max Ceram Shades and Essences. For further information, please refer to the IPS e.max Ceram Instructions for Use.


- IPS e.max Ceram Shades are ready-to-use stains in syringes.
- IPS e.max Ceram Essences are intensively shaded powdered stains, which are mixed with IPS e.max Ceram Glaze and Stain Liquid.

The following steps must be observed:

- Clean the pressed object with a steam cleaner to remove any contaminations and grease residue. Any contamination after cleaning must be prevented.
- **Tip:** For the characterization, place the abutment or the abutment crown on the Ti base using a little IPS e.max Ceram Glaze and Stain Liquid. This allows you to assess the effect of the Ti base on the shade.
- For better wetting of the stains, a small quantity of IPS e.max Ceram Glaze and Stain Liquid may be slightly rubbed into the area that needs to be characterized.
- Mix the pastes or powders with the IPS e.max Ceram Glaze and Stain Liquids (allround or longlife) until the desired consistency is achieved.
- More intensive shades are achieved by several staining procedures and repeated firing, not by applying thicker layers.
- To imitate the incisal area and translucency of the abutment crown in the incisal and occlusal third, IPS e.max Ceram Shade Incisal may be used.
- The cusps and fissures can be individualized using Essences.
- If abutments are fabricated, only the area of the emergence profile is characterized with IPS e.max Ceram Shades and Essences. Do not apply materials to the bonding surface to the crown, as this might compromise the fit and the bond.
- **Important:** Make sure that absolutely no materials are applied to the screw channel and the interface to the Ti base in order not to compromise the fit and bond.
- Conduct the Stain and Characterization firing on a honey-comb firing tray using the stipulated firing parameters. Please observe the special firing parameters.
- Remove restoration from the furnace after completion of the firing cycle (wait for the acoustic signal of the furnace).
- Allow the objects to cool to room temperature in a place protected from draft.
- Do not touch the hot objects with metal tongs.







Apply IPS e.max Ceram Shade Incisal to imitate the incisal area.



Enhance the chroma of the buccal surface.



Individual characterization of the emergence profile with IPS e.max Ceram Essences.



Conduct the Stain and Characterization firing on a honey-comb firing tray.

#### Firing parameters for the Stain and Characterization firing

IPS e.max Ceram on IPS e.max Press Abutment Solutions	B °C	S min	t <sub>r</sub> °C/min	T °C	H min	V <sub>1</sub> °C	V <sub>2</sub> °C	L °C
Stain and Characterization firing	403	6:00	60	770	01:00	450	769	500

Additional Stain and Characterization firing cycles can be conducted with the same firing parameters.

**Note:**

If the layer thickness is less than 2 mm on the entire pressed object, long-term cooling (L) is not required to produce a tension-free condition.

## Glaze firing

Glaze firing is conducted with powder or paste glaze. On abutments, only the emergence profile is glazed. Glaze firing may also take place at a later stage, i.e. when the crown is glazed. On abutment crowns, glaze is applied to the entire outer surface.



The following procedure is recommended:


- Mix the glazing material (IPS e.max Ceram Glaze Paste or Powder) with the IPS e.max Ceram Glaze and Stain Liquids allround or longlife to the desired consistency.
- Apply an even layer of glazing material covering all areas that are to be glazed.
- If required, the fluorescence may be increased by applying a fluorescing glazing material (paste or powder).
- **Important:** Make sure that absolutely no materials are applied to the screw channel and the interface to the Ti base in order not to compromise the fit.
- Make sure that no glaze material is present on the interface of abutments and abutment crowns prior to the firing cycle. If necessary, carefully remove the glaze material.
- Conduct the Glaze firing on a honey-comb firing tray using the stipulated firing parameters. Please observe the special firing parameters.
- Remove restoration from the furnace after completion of the firing cycle (wait for the acoustic signal of the furnace).
- Allow the objects to cool to room temperature in a place protected from draft.
- Do not touch the hot objects with metal tongs.
- If adjustments are required after Glaze firing (e.g. contact points), they may be applied using IPS e.max Ceram Add-on (see page 28).




Apply an even layer of glaze material to the emergence profile of the abutment. Make sure that no glaze material enters the screw channel.



Apply the glazing material evenly on the outer surface of the abutment crown. Make sure that no glaze material enters the screw channel.



Make sure that no glaze material is present on the interface of abutments and abutment crowns prior to the firing cycle. If necessary, carefully remove the glaze material.



Conduct the Glaze firing on a honey-comb firing tray with the corresponding parameters.

**Firing parameters for the Glaze firing**

IPS e.max Ceram on IPS e.max Press Abutment Solutions	B °C	S min	t↗ °C/min	T °C	H min	V <sub>1</sub> °C	V <sub>2</sub> °C	L °C
Glaze firing	403	6:00	60	770	1:00– 2:00	450	769	500

If the gloss is unsatisfactory after the first Glaze firing, further Glaze firing procedures may be conducted using the same firing parameters.

**Note:**

If the layer thickness is less than 2 mm on the entire pressed object, long-term cooling (L) is not required to produce a tension-free condition.



Completely glazed and characterized abutment and abutment crown

**Adjustments with IPS e.max Ceram Add-On**

Use IPS e.max Ceram Add-On Dentin and/or Incisal to make adjustments to the abutment or abutment crown after Glaze firing. Please observe the following procedure for processing:

- Mix IPS e.max Ceram Add-On Dentin or Incisal with IPS e.max Ceram Build-Up Liquid soft or allround and apply on the corresponding areas.
- Fire with the stipulated parameters for the "Add-On after Glaze firing". Observe long-term cooling! If necessary, polish the adjusted areas to a high gloss after firing.



**Firing parameters IPS e.max Ceram Add-On after Glaze firing.**

IPS e.max Ceram on IPS e.max Press	B °C	S min.	t↗ °C/min	T °C	H min	V <sub>1</sub> °C	V <sub>2</sub> °C	L °C
Glaze firing	403	6:00	50	700	01:00	450	699	500

**Optional: Polishing the emergence profile of the abutment**

If no characterizations and no glaze firing of the abutment are required, the emergence profile may be manually polished. Please remember that polishing results in a slight reduction of the emergence profile, which might influence the fit to the gingiva in certain situations.

Please observe the following procedure for polishing:

- Clean the pressed object with a steam cleaner to remove any contaminations.
- Place the pressed object onto the Ti base for processing.
- Overheating of the glass-ceramic must be avoided. Observe the recommendations of the manufacturer of the grinding tool.
- Pre-polish the emergence profile with a diamond-coated rubber polisher. **Note:** The Ti base should not be modified.
- Fine polishing of the emergence profile with a high-gloss rubber polisher.
- High-gloss polishing with brushes and polishing paste.
- Clean the abutment with an ultrasonic or steam cleaner.



Pre-polish the emergence profile with a diamond-coated rubber polisher.



High-gloss polishing with brushes and polishing paste.



Then clean the abutment with ultrasonic ...



... or with steam.



## Crown on hybrid abutment

Preferably, a crown made of IPS e.max Press is seated on an IPS e.max Press hybrid abutment. Depending on your preference, the staining, cut-back or layering technique may be applied for this. For a detailed description of the fabrication, please refer to the IPS e.max Press Instructions for Use.


The following paragraphs describe the working steps which deviate from the procedure applied for the fabrication of restorations on prepared teeth. In this example, an IPS e.max Press crown is fabricated in the cut-back technique.



Screw the Ti base onto the model analog with the corresponding screw. If required, the abutment can be secured on the Ti base by means of Virtual Extra Light Body Fast Set. Seal the screw channel (e.g. with silicone).



Before contouring the crown, apply a spacer up to approx. 1 mm from the cervical crown margin.



Isolate the abutment and then contour the crown with wax. Finally, the restoration is pressed with IPS e.max Press material.



Pressed IPS e.max Press crown with cut-back after divesting and finishing.



Complete the anatomical shape of the reduced areas using IPS e.max Ceram layering materials, such as Incisal and Opal.



Finish the restoration with diamonds and give it a true-to-nature shape and surface structure.



Finally, conduct the Stain and Glaze firing with IPS e.max Ceram Shades, Essences and Glaze.



Abutment and matching crown after Characterization and Glaze firing.

# PS e.max<sup>®</sup> Press Abutment Solutions

## Optional: Clinical Try-in

### Temporarily securing the pressed object on the Ti base

Before the abutment or abutment crown is permanently luted to the Ti base, a clinical try-in may be performed. To facilitate the intraoral handling, the components are temporarily attached to one another with silicone material, e.g. Virtual<sup>®</sup> Extra Light Body Fast Set.

Please observe the following procedure to temporarily secure the components in place:

- Clean the non-pre-treated Ti base and the pressed object (abutment or abutment crown) with steam and blow dry.
- Place the pressed object on the Ti base (which is screwed to the model analog) and mark the relative position of the components. This facilitates the achievement of the correct position when the parts are subsequently temporarily assembled.
- Seal the screw channel with a foam pellet.
- Insert the Virtual cartridge into the dispenser and remove the protective cap.
- Screw on the mixing tip and attach the Oral Tip to the mixing tip.
- Apply Virtual Extra Light Body Fast Set both to the Ti base and directly into the pressed object.
- Insert the Ti base into the pressed object. Observe the relative position of the objects (rotation lock/marking).
- Hold the parts in the correct relative position for 2:30 minutes until the Virtual Extra Light Body Fast Set has set.
- Carefully remove protruding excess material with a suitable instrument, e.g. a scalpel.



Cleaned, non-pre-treated pressed objects (abutment or abutment crown)



Place the abutment or abutment crown onto the Ti base and mark the relative position.



Seal the screw channel with a foam pellet.



Insert the Virtual cartridge into the dispenser, screw on the mixing tip and attach the Oral Tip.



Apply Virtual Extra Light Body Fast Set both to the Ti base ...



... and directly into the pressed object (abutment/abutment crown).



Insert the Ti base into the pressed object. In doing so, observe the relative position of the two components (rotation lock/markings). Hold the components in place for approx. 2:30 minutes until the Virtual Extra Light Body Fast Set material has set.



Carefully remove protruding excess material with a suitable instrument, e.g. a scalpel.



Remove excess Virtual Extra Light Body Fast Set material from the screw channel with an instrument.

## Clinical try-in

### Hybrid abutment

**Important note:** Any intraoral checking of the occlusion/articulation and possible adjustments by grinding may only be performed if the objects have been attached to one another by means of Virtual Extra Light Body Fast Set. During try-in, the Virtual material acts as a buffer, particularly if grinding is necessary, and prevents chipping in the transition area between the hybrid abutment and the crown.

Please observe the following procedure for the clinical try-in:

- Have the clean prepared hybrid abutment (temporarily secured) and the matching clean crown ready at hand.
- Remove the temporary restoration.
- Manually screw in the hybrid abutment with the matching screw.
- Check the geometry of the hybrid abutment (e.g. fit, gingival anaemia) with regard to the gingival margin.
- If required, seal the screw channel on the hybrid abutment with a foam pellet.
- **Tip:** Isolate the inner aspect of the crown with glycerin gel, e.g. Try-in paste, Liquid Strip
- **Place the crown intraorally onto the hybrid abutment to check and adjust the proximal contacts, if necessary.**
- **Note: No occlusal functional checks must be performed at this stage.**
- **For the functional check, the crown has to be secured on the hybrid abutment with Virtual Extra Light Body Fast Set. Do not use Try-in paste for this purpose, as this material is not sufficiently resistant against the compressive forces.**
- Insert the Virtual cartridge into the dispenser and remove the protective cap.
- Screw on the mixing tip and attach the Oral Tip to the mixing tip.
- Apply Virtual Extra Light Body Fast Set to the inner aspect of the crown.
- Use your finger to press the crown to the hybrid abutment until the final position has been achieved. Hold the crown in the final position until the Virtual material has set.
- Remove excess Virtual material.
- Check the occlusion/articulation and make required adjustments with suitable grinding instruments (see separate IPS e.max Recommended grinding instruments for ceramics – use in the dental practice). If adjustments have been made by grinding, conduct another polishing cycle or glaze firing.
- Carefully remove the crown from the hybrid abutment and the hybrid abutment (including Ti base).
- Insert the temporary restoration.



Manually screw in the hybrid abutment with the matching screw. Check the geometry of the hybrid abutment (e.g. fit, gingival anaemia) with regard to the gingival margin.



If required, seal the screw channel of the hybrid abutment with a foam pellet.






Tip: Isolate the inner aspect of the crown with glycerine gel.



Place the crown intraorally onto the hybrid abutment to check and adjust the proximal contacts, if necessary. **Note: No occlusal functional checks must be performed at this stage.**



Apply Virtual Extra Light Body Fast Set to the inner aspect of the crown.



Use your finger to press the crown to the hybrid abutment until the final position has been achieved. Hold the crown in the final position until the Virtual material has set.



Remove excess Virtual material.



Check the occlusion/articulation and use suitable grinding instruments to make possibly required adjustments.



Carefully remove the crown from the hybrid abutment and remove the Virtual Extra Light Body Fast Set material.



Unscrew the hybrid abutment.

### Hybrid abutment crown

Please observe the following procedure for the clinical try-in:

- Have the cleaned hybrid abutment crown (temporarily secured with Virtual Extra Light Body Fast Set) at hand.
- Remove the temporary restoration.
- **Place the hybrid abutment crown intraorally onto the implant to check and possibly adjust the proximal contacts. Note: No occlusal functional checks must be performed at this stage.**
- Manually screw in the hybrid abutment crown with the matching screw.
- Check the geometry of the hybrid abutment crown (e.g. fit, gingival anaemia) with regard to the gingiva.
- Check the occlusion/articulation and make possibly required adjustments with suitable grinding instruments (see separate IPS e.max Recommended grinding instruments for ceramics – use in the dental practice). If adjustments have been made by grinding, conduct another polishing cycle or glaze firing.
- Carefully remove the hybrid abutment crown (including Ti base).
- Rinse the implant site e.g. with Cervitec Liquid (antibacterial mouth wash with chlorhexidine) to clean and disinfect it.
- Insert the temporary restoration.



Place the hybrid abutment crown intraorally onto the implant to check and possibly adjust the proximal contacts. **Note: No occlusal functional checks must be performed at this stage.**



Manually screw in the hybrid abutment crown with the matching screw.



Check the geometry of the hybrid abutment crown (e.g. fit, gingival anaemia) with regard to the gingiva.



Check the occlusion/articulation and use suitable grinding instruments to make possibly required adjustments.



Carefully remove the hybrid abutment crown (including Ti base).

# IPS e.max<sup>®</sup> Press Abutment Solutions

## Permanent Cementation

Careful preparation of the bonding surface is a prerequisite for an optimum adhesive cementation of the Ti base and the pressed object. The following paragraphs outline the required procedures. The procedure is the same for hybrid abutments and hybrid abutment crowns.

	IPS e.max Press Abutment Solutions	
	Abutment, abutment crown	Ti base
<b>Blasting</b>	–	The bonding area with Al <sub>2</sub> O <sub>3</sub> at low pressure
<b>Etching</b>	The bonding area with IPS <sup>®</sup> Ceramic Etching Gel for 20 s	–
<b>Conditioning/silanating</b>	The bonding area with Monobond <sup>®</sup> Plus for 60 s	
<b>Adhesive cementation</b>	Multilink <sup>®</sup> Implant MO 0	
<b>Covering the cementation joint</b>	Glycerine gel, e.g. Liquid Strip	
<b>Curing</b>	7-minute curing (optionally in a light-curing device)	
<b>Polishing the cementation joint</b>	Customary polishers for ceramic/resin materials	

All the materials that are required for permanent cementation and for clinical try-in are contained both in the **IPS e.max Press Abutment Solutions Basic Kit A–D\*** and the **IPS e.max Abutment Solutions Cem Kit\***.



\* The range of available products may vary from country to country.



## Pre-treatment of the Ti base

To prepare the Ti base for cementation with the pressed object, please observe the following procedure:


- Observe processing instructions of the manufacturer of the Ti base.
- Clean the Ti base in an ultrasonic cleaner and blow dry or use a steam cleaner.
- Screw Ti base onto a model analog.
- Place the pressed object on the Ti base and mark the relative position of the components. This facilitates locating the correct position when the parts are assembled at a later stage.
- The emergence profile of the Ti base must not be blasted or modified in any way. To protect the emergence profile, hard modelling wax is applied, as this material can be easily removed later on.
- Seal the screw channel with wax.
- Carefully blast the bonding area with  $\text{Al}_2\text{O}_3$  (50–100  $\mu\text{m}$ ) at low pressure until an even mat surface has been achieved.
- Clean with an instrument and steam cleaner. Make sure that any wax residue is carefully removed.
- After cleaning, any contamination of the bonding surface must be prevented, since contaminations negatively influence the bond.
- Apply Monobond Plus on the cleaned bonding surface and allow to react for 60 seconds. After the reaction time, dry the remaining residue with water- and oil-free air.
- Seal the screw channel with a foam pellet or wax. Make sure that the bonding surface is not contaminated.



Screw Ti base onto a model analog. Mark the relative position to the pressed object.



Apply wax to protect the emergence profile. In addition, seal also the screw channel with wax.



Carefully blast the bonding area with  $\text{Al}_2\text{O}_3$  (50–100  $\mu\text{m}$ ) at low pressure until an even mat surface has been achieved.



Clean with an instrument and steam cleaner. Make sure that any wax residue is carefully removed.



Apply Monobond Plus on the cleaned bonding surface and allow to react for 60 seconds. After the reaction time, dry the remaining residue with water- and oil-free air.



Seal the screw channel with a foam pellet or wax. Make sure that the bonding surface is not contaminated.

## Preparing the pressed object

To prepare the pressed object for cementation onto the Ti base, please observe the following procedure:

- Do **not** blast the IPS e.max Press object to prepare it for cementation.
- Clean the pressed object in an ultrasonic cleaner and blow dry or use a steam cleaner.
- After cleaning, any contamination of the bonding surface must be prevented, since contaminations negatively influence the bond.
- To protect the outer surfaces or the glazed areas, wax may be applied.
- Etch the bonding surface with 5% hydrofluoric acid gel (IPS Ceramic Etching Gel) for 20 seconds.
- Subsequently, thoroughly rinse the bonding surface under running water and dry with oil-free air.
- Apply Monobond Plus on the cleaned bonding surface and allow to react for 60 seconds. After the reaction time, dry the remaining residue with water and oil-free air.



Do **not** blast IPS e.max Press object



Etch with IPS Ceramic Etching Gel for 20 seconds.



Allow Monobond Plus to react for 60 seconds and dry with air.

## Cementation with Multilink® Implant

For an optimum bond between the IPS e.max Press object and the Ti base, use the self-curing luting composite Multilink® Implant with light-curing option. Read the respective Instructions for Use for more detailed information.

Please observe the following procedure for cementation:

- Keep the cleaned and conditioned components that are to be luted (pressed objects, Ti base) at hand.
- **The subsequent cementation procedure must be carried out quickly and without interruption. The working time of Multilink Implant is 90 (± 15) seconds at 23 °C (± 1 °C) or 73 °F (± 1.8 °F).**
- As a general rule, attach a new mixing tip to the Multilink Implant syringe prior to each use.
- Apply Multilink Implant directly from the mixing tip in a thin layer on the bonding surface of the **Ti base and the bonding surface of the pressed object**.
- Leave the mixing tip attached to the Multilink Implant syringe until the next application. As the material will cure in the mixing tip, it will serve as a seal.
- Position the pressed object above the Ti base in such a way that the position markings are aligned.
- Use even and low pressure to join the parts and check the correct relative position of the components when they are in their final position (transition Ti base/pressed object).
- Subsequently, tightly press the components together for 5 seconds.
- Carefully remove excess in the screw cavity, e.g. with Microbrush or brush, using a rotary movement.
- Remove excess at the transition to the Ti base carefully in its ductile state, e.g. with a foam pellet, while applying slight pressure to hold the components in place.
- Apply a glycerine gel (e.g. Liquid Strip) to the cementation joint to prevent the formation of a inhibition layer.
- After that, completely polymerize the luting composite for 7 minutes in a light-polymerization device.
- **Important: Do not move the objects until the Multilink Implant material has completely cured and hold them in place without allowing for any motion, e.g. with diamond-coated tweezers.**
- After completion of the polymerization, rinse off the glycerine gel with water.
- **Smooth out and polish the cementation joint with rubber polishers.**
- If there are any cement residues in the screw channel, remove them using suitable rotary instruments.
- Clean with a steam cleaner.




Keep the cleaned and conditioned components that are to be luted at hand.



Attach a new mixing tip to the Multilink Implant syringe prior to each use.




Apply Multilink Implant directly from the mixing tip in a thin layer on the bonding surface of the Ti base.



Apply Multilink Implant directly from the mixing tip in a thin layer on the bonding surface of the pressed object.



Position the pressed object above the Ti base in such a way that the position markings are aligned.



Join the components using even and light pressure. Subsequently, tightly press the components together for 5 seconds.



Carefully remove excess in the screw cavity, e.g. with a Microbrush or a brush, using a rotary movement.



Remove excess carefully in its ductile state, e.g. with a foam pellet, while applying slight pressure to hold the components in place.



Apply a glycerine gel (e.g. Liquid Strip) to the cementation joint to prevent the formation of an inhibition layer.



Polymerize the luting composite for 7 minutes (optionally in a light-curing unit). **Important: Do not move the objects until the material has completely cured and hold them in place without allowing for any motion.**



After completion of the polymerization, rinse off the glycerine gel with water.



Smooth out and polish the cementation joint with rubber polishers.



If there is any cement residue in the screw channel, remove using suitable rotary instruments. Do not damage the Ti base.



Completed hybrid abutment and hybrid abutment crown after cementation.



# IPS e.max<sup>®</sup> Press Abutment Solutions

## Seating and Aftercare

### Sterilization

We recommend sterilizing hybrid abutments or hybrid abutment crowns before seating them in the patient's mouth.

- The sterilization time is 15 minutes at 121 °C / 250 °F.
- Only devices which comply with the standards EN 13060 and EN 285 should be used for sterilization. The sterilization processes are validated according to EN ISO-17664:2004.



### Intraoral preparation

Please observe the following procedure to prepare for the permanent cementation of the implant-supported restoration:

- Remove the temporary restoration.
- Clean the implant site.
- Check the periimplant tissue (emergence profile).

### Seating the hybrid abutment and crown

#### Preparing/conditioning the hybrid abutment and crown

Conditioning of the ceramic surface, i.e. the bonding surface, in preparation for cementation is critical for generating a sound bond between the cementation material and the all-ceramic material.

The following steps must be observed:

- Do not blast the IPS e.max Press hybrid abutment or the IPS e.max Press crown with Al<sub>2</sub>O<sub>3</sub> or glass polishing beads prior to seating.
- Ideally, conduct the clinical try-in before etching in order not to contaminate the bonding surface.
- Thoroughly clean the hybrid abutment and the crown with water and blow dry.
- Etch the bonding surfaces with 5% hydrofluoric acid gel (IPS Ceramic Etching Gel) for 20 seconds. Make sure that no etching gel comes into contact with the emergence profile or the outer side of the crown. **Important: Do not use the IPS Ceramic Etching Gel intraorally.**
- Thoroughly rinse off the etching gel with water and dry with oil- and water-free air.
- If an adhesive or self-adhesive cementation protocol is used, apply Monobond Plus to the bonding surfaces, allow to react for 60 seconds and then dry with oil- and water-free air.



Do not blast IPS e.max Press objects.



Etch the bonding surfaces with IPS Ceramic Etching Gel for 20 seconds.



Apply Monobond Plus to the bonding areas, allow to react for 60 s and blow dry.

### Seating the hybrid abutment and crown

For the permanent seating of the hybrid abutment and the crown, please observe the following working steps. Please also observe the Instructions for Use of the selected luting material.


- Do not use phenolic mouth washes, as such products negatively influence the bond between the ceramic and the composite.
- Insert the hybrid abutment intraorally into the implant.
- Manually screw in the matching implant screw.
- Tighten the implant screw with a torque wrench (observe the instructions of the manufacturer).
- Insert a cotton or foam pellet into the screw channel.
- Seal the screw channel with a temporary composite (e.g. Telio® CS Inlay). This serves to ensure access to the screw at a later stage.
- Check the bonding area for contamination/moisture and clean or dry with an air syringe, if necessary.
- Apply the luting material, e.g. SpeedCEM®, into the conditioned crown.
- Place the crown onto the hybrid abutment and secure in place in the final position.
- Conduct the pre-polymerization using the four-quarter technique.
- Remove excess luting material.
- Cover the cementation joint with glycerine gel (e.g. Liquid Strip).
- Polymerize with an LED curing light (e.g. bluephase®).
- Rinse off the glycerine gel with water.
- Check the occlusion and articulation and make adjustments, if necessary. If adjustments are made to the restoration by grinding, these areas must subsequently be polished to a high gloss, e.g. using OptraFine.
- Polish restoration margins and the cementation joint with silicone polishers (e.g. Astropol®, OptraFine).
- Apply Cervitec Plus in the area of the gingival margin.



Insert the hybrid abutment intraorally into the implant.



Manually screw in the matching implant screw.



Tighten the implant screw with a torque wrench (observe the instructions of the manufacturer).




Seal the screw channel, for instance with a cotton or foam pellet and a temporary composite material.



Apply the luting material, e.g. SpeedCEM, into the conditioned crown.



Place the crown onto the hybrid abutment and secure in place.



Conduct the pre-polymerization using the four-quarter technique.



Remove excess luting material.



Cover the restoration margin with glycerine gel (e.g. Liquid Strip).



Polymerize with an LED curing light (e.g. bluephase).



Rinse off the glycerine gel with water.



Check the occlusion and articulation and make adjustments, if necessary.





Polish restoration margins and the cementation joint with polishers (e.g. Astropol, OptraFine).



Completed IPS e.max Press hybrid abutment and crown.



## Seating the hybrid abutment crown.

### Preparing/conditioning the hybrid abutment crown

Please observe the following notes to prepare for the intraoral sealing of the screw channel:

- As a general rule, do **not** blast IPS e.max Press hybrid abutment crowns with  $\text{Al}_2\text{O}_3$  or glass polishing beads.
- Thoroughly clean the the hybrid abutment crown with water and blow dry.
- Etch the screw channel from the occlusal side with 5% hydrofluoric acid gel (IPS Ceramic Etching Gel) for 20 seconds. Make sure that no etching gel comes into contact with the occlusal surface. **Important: Do not use the IPS Ceramic Etching Gel intraorally.**
- **Thoroughly rinse off the etching gel with water and dry with oil- and water-free air.**
- Apply Monobond Plus to the etched and cleaned surface in the screw channel, allow to react for 60 seconds and then blow dry with oil- and water-free air.



Do **not** blast IPS e.max Press object



Etch with IPS Ceramic Etching Gel for 20 seconds.



Apply Monobond Plus and allow to react for 60 seconds and blow dry with air.

### Seating the hybrid abutment crown.

For the permanent seating of the hybrid abutment crown, please observe the following working steps.

- Do not use phenolic mouth washes, as such products negatively influence the bond between the ceramic and the composite.
- Insert the hybrid abutment crown intraorally into the implant.
- Manually screw in the matching implant screw.
- Tighten the implant screw with a torque wrench (observe the instructions of the manufacturer).
- Check the screw channel for contamination/moisture and clean with Total Etch (phosphoric acid gel), if necessary.
- Insert a cotton or foam pellet into the screw channel.
- Apply the bonding agent.
- Seal the screw channel with a composite material (e.g. Tetric EvoCeram) in the appropriate shade.
- Polymerize with an LED curing light (e.g. bluephase).
- Check the occlusion/articulation after polymerization and correct possible rough spots with suitable fine-grain diamonds.
- Polish to a high gloss with silicone polishers (e.g. OptraFine).



Insert the hybrid abutment crown intraorally into the implant.



Manually screw in the matching implant screw.



Tighten the implant screw with a torque wrench (observe the instructions of the manufacturer).



Seal the screw channel with a composite material (e.g. Tetric EvoCeram) in the appropriate shade.



Polymerize with an LED curing light (e.g. bluephase).



After polymerization, check the occlusion/articulation and correct possible rough spots with suitable finishers (e.g. Astropol F) or fine diamonds.



Polish to a high gloss using silicone polishers (e.g. Astropol P, Astropol HP or Astrobrush).



Completed IPS e.max Press hybrid abutment crown

## Care Notes – Implant Care

Implant Care comprises a range of coordinated products for the professional care of patients during the various phases of implant treatment and lifelong aftercare. Products for professional tooth cleaning and bacterial control contribute to the long-term quality assurance of implant-supported restorations. Structural elements, peri-implant tissue, natural teeth, dental restorations, gingiva and the mucosa are treated in an optimum way with regard to function and esthetics.



# IPS e.max<sup>®</sup> Press Abutment Solutions

## General Information

### Frequently Asked Questions

**In addition to the desired tooth shade, why should the root shade also be defined/determined upon shade determination?**

*IPS e.max Press Abutment Solutions allow you to fabricate restorations with a lifelike appearance both in the visible area and the area below the gingiva (root). By defining the root shade, a highly esthetic outcome can be achieved especially receding gingiva.*

**Is it possible to fabricate an abutment or an abutment crown with IPS e.max Press (LS<sub>2</sub>) without using a Ti base?**

*No! For this indication, IPS e.max Press needs the support provided by the Ti base. In addition, the Ti base allows an optimum (industrially fabricated) fit to the implant to be achieved.*

**Is it possible to use any commercially available Ti base in conjunction with IPS e.max Press Abutment Solutions?**

*When selecting a suitable Ti base, the requirements in terms of minimum dimensions (height, shoulder width, no undercuts) must be taken into account. In addition, the Ti base must be equipped with a rotation lock which does not entail a reduction of the ceramic layer thickness.*

**Is it permissible to modify the selected Ti base?**

*The instructions of the manufacturer with regard to modifying the Ti base must be observed. Prior to permanent cementation, the bonding surface of the Ti base must be blasted with Al<sub>2</sub>O<sub>3</sub>.*

**Is a hybrid abutment crown indicated in the anterior region?**

*This indication depends on the position and inclination of the implant. If the opening of the screw channel is located on the lingual/palatinal surface, a hybrid abutment crown may be fabricated in the anterior region.*

**Is it possible to use an IPS e.max Press hybrid abutment as an abutment for a bridge restoration?**

*No. Only single-tooth restorations may be fabricated.*

**What do I need to take into consideration when designing a hybrid abutment or hybrid abutment crown in order to fabricate a durable restoration?**

*The stipulated minimum and maximum layer thicknesses for IPS e.max Press need to be observed. In addition, the ratio between the height of the Ti base and the height of the entire restoration must be observed.*

**What do I need to take into consideration when attaching sprues and investing the wax-up?**

*The screw channel of the sprued wax-up must be parallel to the outer wall of the investment ring. As a result, the investment material can be filled evenly and in a controlled manner. In addition, the risk that the flowing ceramic material could break off the investment material in the screw channel is reduced. The objects could be placed in a tilted position on the investment ring base, but this may lead to difficulties during investing (e.g. bubbles in the screw channel).*

*If the screw channel is very long, the investment material in the screw channel can be stabilized with a pin (e.g. high-grade steel, ZrO<sub>2</sub>) during investing. For this purpose, pour investment material into the investment ring up to the restoration margin, insert the pin into the screw channel and fill up the investment ring up to the marking without vibrating it.*

**When is the sprue of the pressed object separated?**

*We recommend that the pressed objects should be fitted to the Ti base first, as this facilitates the handling. Subsequently, the pressed objects are separated from the sprue.*

**How should the emergence profile of the hybrid abutment be finished?**

Preferably, a Characterization and Glaze firing is conducted on the emergence profile prior to the cementation procedure. In this way, you may adjust the esthetic appearance of the abutment to the clinical situation ("root shade"). If no characterization is required, the emergence profile may be polished to a high gloss with polishers as well as brushes and paste.

**Is it possible to use IPS e.max Ceram Glaze Spray to glaze the abutment or the abutment crown?**

We do not recommend using the Glaze Spray in this indication, as the bonding surface and the screw channel may be contaminated with glaze material.

**An optional clinical try-in may be performed. How are the objects prepared for this?**

The Ti base and the pressed abutment or abutment crown are temporarily joined in the laboratory by means of a silicone material, e.g. Virtual Extra Light Body Fast Set. This facilitates the intraoral handling.

**What must be observed for the clinical try-in of a crown on a hybrid abutment?**

To check the occlusion/articulation and to make possible adjustments, the crown must be temporarily secured on the hybrid abutment with a silicone material, e.g. Virtual Extra Light Body Fast Set. The silicone material acts as a buffer and prevents chipping in the marginal area of the crown. Try-in pastes or Vaseline must not be used for functional checks.

**What material is used to permanently cement the abutment or the abutment crown made of IPS e.max Press with the Ti base?**

Only Multilink Implant is to be used for the permanent cementation. Other luting materials have not been tested for this purpose.

**How is the Ti base prepared for the permanent cementation with Multilink Implant?**

Carefully blast the bonding area with  $Al_2O_3$  at low pressure until an even mat surface has been achieved. After cleaning, the area is conditioned with Monobond Plus.

**How is the screw channel of a hybrid abutment crown sealed intraorally?**

After the hybrid abutment crown has been screwed into place and the screw has been tightened with a torque wrench, the screw channel is sealed with a composite restorative material.



# Material Selection Table

## Hybrid abutment and separate crown

The material is selected on the basis of the desired tooth shade (Bleach BL or A-D). Depending on the Ti base selected and the design of the hybrid abutment or crown, characterizations with IPS e.max Ceram Shades and Essences may be necessary to achieve the desired shade. The desired tooth shade is achieved after the restoration has been seated and is made up of the shade of the hybrid abutment and the shade of the crown that is cemented onto it. The ingot recommendations for the hybrid abutment have been selected in such a way that the desired tooth shade is achieved in combination with the crown. In the "cervical area", it may be necessary to characterize the hybrid abutment according to the clinical situation.

Desired tooth shade: Bleach BL and A-D Shade Guide																				
	BL1	BL2	BL3	BL4	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Ti base																				
Multilink Implant MO 0																				
Cementation in the laboratory																				
MO* for the hybrid abutment	-	MO 0	MO 0	MO 0	MO 0	MO 0	MO 0	MO 0	MO 0	MO 0	MO 0	MO 0	MO 0	MO 0	MO 0	MO 0	MO 0	MO 0	MO 0	MO 0
HO	HO 0	HO 1	HO 1	HO 1	HO 1	HO 1	HO 1	HO 1	HO 1	HO 1	HO 1	HO 1	HO 1	HO 1	HO 1	HO 1	HO 1	HO 1	HO 1	HO 2
Cementation (intraoral)	adhesive, self-adhesive or conventional cementation e.g. Multilink Implant, SpeedCEM, etc.																			
IPS e.max Press crown	LT BL1	LT BL2	LT BL3	LT BL4	LT A1	LT A2	LT A3	LT A3.5	LT A4	LT B1	LT B2	LT B3	LT B4	LT C1	LT C2	LT C3	LT C4	LT D2	LT D3	LT D4

\* Given the high opacity, HO ingots are better suited to cover the (grey) titanium bases. These ingots may also be used for thinner hybrid abutment designs.

## Hybrid abutment crown



The material is selected on the basis of the desired tooth shade (Bleach BL or A-D). Depending on the Ti base selected and the design of the hybrid abutment crown, characterization with IPS e.max Ceram Shades and Essences may be necessary to achieve the desired shade.



Desired tooth shade: Bleach BL and A-D Shade Guide																				
	BL1	BL2	BL3	BL4	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Ti base																				
Multilink Implant MO 0																				
Cementation in the laboratory																				
Abutment crown IPS e.max Press	LT BL1	LT BL2	LT BL3	LT BL4	LT A1	LT A2	LT A3	LT A3.5	LT A4	LT B1	LT B2	LT B3	LT B4	LT C1	LT C2	LT C3	LT C4	LT D2	LT D3	LT D4

## Press and Firing Parameters

### Press parameters for IPS e.max Press

The press furnace, investment ring size and the selected IPS e.max Press ingot must be considered:

Press furnace	IPS e.max Press ingots	IPS Investment Ring System	B °C	t <sub>r</sub> °C/min	T °C	H min	V <sub>1</sub> °C	V <sub>2</sub> °C	
EP 500		100 g	700	60	925	15	500	925	Program 11-20 Software 2.9
		200 g	700	60	930	25	500	930	Program 11-20 Software 2.9
		100 g	700	60	920	15	500	920	Program 11-20 Software 2.9
		200 g	700	60	925	25	500	925	Program 11-20 Software 2.9

Press furnace	IPS e.max Press ingots	IPS Investment Ring System	B °C	t <sub>r</sub> °C/min	T °C	H min	A
EP 600 Combi		100 g	700	60	915	15	300 µm/min
		200 g	700	60	920	25	300 µm/min
		100 g	700	60	910	15	300 µm/min
		200 g	700	60	915	25	300 µm/min

#### Programat EP 3000



Select the press program in accordance with the selected ingot to be pressed and the investment ring size used.



The press parameters for HO, MO, LT and HT are integrated starting with software V 6.1.

#### Programat EP 5000



Select the press program in accordance with the selected ingot to be pressed and the investment ring size used.



The press parameters for HO, MO, LT and HT are integrated starting with software V 6.1.

- The firing parameters listed represent standard values and apply to the Ivoclar Vivadent furnaces P300, P500, P700, EP 3000 and EP 5000. The temperatures indicated also apply to furnaces of older generations. However, the temperature in the firing chamber may deviate by approx. ± 10 °C/18 °F, depending on the age of the heating muffle.
- If a non-Ivoclar Vivadent furnace is used, temperature corrections may be necessary.
- Regional differences in the power supply or the operation of several electronic devices by means of the same circuit may also render adjustments of the temperatures necessary.



### Firing parameters for IPS e.max Press Abutment Solutions

- Use a honey-comb tray and the corresponding pins for firing.
- Do not use ceramic pins.
- The stipulated parameters are coordinated with Ivoclar Vivadent furnaces (tolerance range +/- 10 °C/18 °F).
- If furnaces other than those from Ivoclar Vivadent are used, temperature adjustments may be necessary.
- Remove IPS e.max Press objects from the furnace after completion of the firing cycle (wait for the acoustic signal of the furnace).
- Allow the objects to cool to room temperature in a place protected from draft.
- Do not touch the hot objects with metal tongs.
- Do not blast or quench the objects.

### Note regarding the firing parameters for IPS e.max Press Abutment Solutions:

Given the geometry of hybrid abutments and hybrid abutment crowns, the layer thicknesses of the pressed object may vary considerably. When the objects cool after the firing cycle, the different cooling speeds in the areas with different thicknesses may result in internal tensions. In the worst case, these tensions may result in fractures in the pressed ceramic objects. By using slow cooling (long-term cooling (L)), these tensile stresses can be minimized.

### Note:

Please read the notes of the manufacturer of the ceramic furnace for more details on the programming of the long-term cooling (L). If the layer thickness is less than 2 mm on the entire pressed object, long-term cooling (L) is not required.



IPS e.max Ceram on IPS e.max Press Abutment Solutions	B °C	S min	t <sup>↗</sup> °C/min	T °C	H min	V <sub>1</sub> °C	V <sub>2</sub> °C	L °C
Stain and Characterization firing	403	6:00	60	770	01:00	450	769	500
Glaze firing	403	6:00	60	770	1:00 – 2:00	450	769	500
Add-On after Glaze firing	403	6:00	50	700	01:00	450	699	500

## Clinical cases (R. Watzke, Liechtenstein)

IPS e.max Press hybrid abutment / IPS e.max Press crown (area 36), IPS e.max Press crown (area 37).



Clinical situation after implantation and gingiva shaping



IPS e.max Press hybrid abutment (36) and IPS e.max Press crown, IPS e.max Press single crown



IPS e.max Press hybrid abutment, screwed in, cemented IPS e.max Press crown



Final image, buccal view: IPS e.max Press crown cemented onto an IPS e.max Press hybrid abutment



Final image, occlusal view: IPS e.max Press crown cemented onto an IPS e.max Press hybrid abutment

IPS e.max Press hybrid abutment crown (area 35)



Clinical situation after implantation and gingiva shaping



IPS e.max Press hybrid abutment crown



Screwing in the hybrid abutment crown



Sealing the screw channel with composite material (e.g. Tetric EvoCeram)



Final image of an IPS e.max Press hybrid abutment crown

# Ivoclar Vivadent – worldwide

## **Ivoclar Vivadent AG**

Bendererstrasse 2  
FL-9494 Schaan  
Liechtenstein  
Tel. +423 235 35 35  
Fax +423 235 33 60  
www.ivoclarvivadent.com

## **Ivoclar Vivadent Pty. Ltd.**

1 – 5 Overseas Drive  
P.O. Box 367  
Noble Park, Vic. 3174  
Australia  
Tel. +61 3 979 595 99  
Fax +61 3 979 596 45  
www.ivoclarvivadent.com.au

## **Ivoclar Vivadent Ltda.**

Alameda Caiapós, 723  
Centro Empresarial Tamboré  
CEP 06460-110 Barueri – SP  
Brazil  
Tel. +55 11 2424 7400  
Fax +55 11 3466 0840  
www.ivoclarvivadent.com.br

## **Ivoclar Vivadent Inc.**

2785 Skymark Avenue, Unit 1  
Mississauga  
Ontario L4W 4Y3  
Canada  
Tel. +1 905 238 5700  
Fax +1 905 238 5711  
www.ivoclarvivadent.com

## **Ivoclar Vivadent Marketing Ltd.**

Rm 603 Kuen Yang  
International Business Plaza  
No. 798 Zhao Jia Bang Road  
Shanghai 200030  
China  
Tel. +86 21 5456 0776  
Fax +86 21 6445 1561  
www.ivoclarvivadent.com

## **Ivoclar Vivadent Marketing Ltd.**

Calle 134 No. 7-B-83, Of. 520  
Bogotá  
Colombia  
Tel. +57 1 627 33 99  
Fax +57 1 633 16 63  
www.ivoclarvivadent.co

## **Ivoclar Vivadent SAS**

B.P. 118  
F-74410 Saint-Jorioz  
France  
Tel. +33 450 88 64 00  
Fax +33 450 68 91 52  
www.ivoclarvivadent.fr

## **Ivoclar Vivadent GmbH**

Dr. Adolf-Schneider-Str. 2  
D-73479 Ellwangen, Jagst  
Germany  
Tel. +49 (0) 79 61 / 8 89-0  
Fax +49 (0) 79 61 / 63 26  
www.ivoclarvivadent.de

## **Ivoclar Vivadent Marketing (India) Pvt. Ltd.**

503/504 Raheja Plaza  
15 B Shah Industrial Estate  
Veera Desai Road, Andheri (West)  
Mumbai, 400 053  
India  
Tel. +91 (22) 2673 0302  
Fax +91 (22) 2673 0301  
www.ivoclar-vivadent.in

## **Ivoclar Vivadent s.r.l.**

Via Isonzo 67/69  
40033 Casalecchio di Reno (BO)  
Italy  
Tel. +39 051 611 35 55  
Fax +39 051 611 35 65  
www.ivoclarvivadent.it

## **Ivoclar Vivadent K.K.**

1-28-24-4F Hongo  
Bunkyo-ku  
Tokyo 113-0033  
Japan  
Tel. +81 3 6903 3535  
Fax +81 3 5844 3657  
www.ivoclarvivadent.jp

## **Ivoclar Vivadent Ltd.**

12F W-Tower, 1303-37  
Seocho-dong, Seocho-gu,  
Seoul 137-855  
Republic of Korea  
Tel. +82 (2) 536 0714  
Fax +82 (2) 596 0155  
www.ivoclarvivadent.co.kr

## **Ivoclar Vivadent S.A. de C.V.**

Av. Insurgentes Sur No. 863.  
Piso 14, Col. Napoles  
03810 México, D.F.  
México  
Tel. +52 (55) 50 62 10 00  
Fax +52 (55) 50 62 10 29  
www.ivoclarvivadent.com.mx

## **Ivoclar Vivadent Ltd.**

12 Omega St, Albany  
PO Box 5243 Wellesley St  
Auckland, New Zealand  
Tel. +64 9 914 9999  
Fax +64 9 814 9990  
www.ivoclarvivadent.co.nz

## **Ivoclar Vivadent Polska Sp. z o.o.**

Al. Jana Pawla II 78  
00-175 Warszawa  
Poland  
Tel. +48 22 635 54 96  
Fax +48 22 635 54 69  
www.ivoclarvivadent.pl

## **Ivoclar Vivadent Marketing Ltd.**

Derbenevskaja Nabereshnaya 11, Geb. W  
115114 Moscow  
Russia  
Tel. +7 495 913 66 19  
Fax +7 495 913 66 15  
www.ivoclarvivadent.ru

## **Ivoclar Vivadent Marketing Ltd.**

Qlaya Main St.  
Siricon Building No.14, 2<sup>nd</sup> Floor  
Office No. 204  
P.O. Box 300146  
Riyadh 11372  
Saudi Arabia  
Tel. +966 1 293 83 45  
Fax +966 1 293 83 44  
www.ivoclarvivadent.com

## **Ivoclar Vivadent Pte. Ltd.**

171 Chin Swee Road  
#02-01 San Centre  
Singapore 169877  
Tel. +65 6535 6775  
Fax +65 6535 4991  
www.ivoclarvivadent.com

## **Ivoclar Vivadent S.L.U.**

c/ Emilio Muñoz Nº 15  
Entrada c/ Albarracín  
E-28037 Madrid  
Spain  
Tel. + 34 91 375 78 20  
Fax + 34 91 375 78 38  
www.ivoclarvivadent.es

## **Ivoclar Vivadent AB**

Dalvägen 14  
S-169 56 Solna  
Sweden  
Tel. +46 (0) 8 514 93 930  
Fax +46 (0) 8 514 93 940  
www.ivoclarvivadent.se

## **Ivoclar Vivadent Liaison Office**

: Tesvikiye Mahallesi  
Sakayik Sokak  
Nisantas' Plaza No:38/2  
Kat:5 Daire:24  
34021 Sisli – Istanbul  
Turkey  
Tel. +90 212 343 08 02  
Fax +90 212 343 08 42  
www.ivoclarvivadent.com

## **Ivoclar Vivadent Limited**

Ground Floor Compass Building  
Feldspar Close  
Warrens Business Park  
Enderby  
Leicester LE19 4SE  
United Kingdom  
Tel. +44 116 284 78 80  
Fax +44 116 284 78 81  
www.ivoclarvivadent.co.uk

## **Ivoclar Vivadent, Inc.**

175 Pineview Drive  
Amherst, N.Y. 14228  
USA  
Tel. +1 800 533 6825  
Fax +1 716 691 2285  
www.ivoclarvivadent.us

Date information prepared: 11/2011

Some products and/or indications may not be regulatory cleared/released in all markets. Please contact the local Ivoclar Vivadent sales office for the current national status.

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

Printed in Liechtenstein  
© Ivoclar Vivadent AG, Schaan / Liechtenstein  
641315/1111/e/BVD

Traditional 510(k)  
Straumann® Variobase™ Abutments

Straumann USA, LLC  
January 23, 2014

  
ivoclar  
vivadent®  
technical

A16-57

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 17

---

**Appendix 17 – IPS e.max® CAD Abutment Solutions Instructions  
for Use**

# ® e.max CAD

Abutment Solutions

e.max  
IPS



Instructions for Use

CE 0123

Traditional 510(k)  
Straumann® Variobase™ Abutments

Straumann USA, LLC  
January 23, 2014

ivoclar  
vivadent  
technical

A17-2

# Table of Contents

<b>3</b>	<b>IPS e.max System</b>
<b>4</b>	<b>IPS e.max CAD</b>
<b>5</b>	<b>Product Information</b> Description of IPS e.max CAD Abutment Solutions Material Scientific data CAD/CAM partners
<b>9</b>	<b>Fabrication of Hybrid Abutments and Hybrid Abutment Crowns</b> Treatment/Fabrication process Shade – tooth shade and abutment shade Preparation for the CAD/CAM process Layer thicknesses of the ceramic components Block selection CAD/CAM processing Finishing
<b>17</b>	<b>Optional: Clinical Try-In</b> Temporarily securing the ceramic structure on the Ti base Clinical try-in
<b>22</b>	<b>Completing the IPS e.max CAD Ceramic Structure</b> Polishing technique Staining technique on the "blue restoration" Staining technique on the "tooth-coloured restoration"
<b>44</b>	<b>Crowns on IPS e.max CAD Hybrid Abutments</b>
<b>46</b>	<b>Permanent Cementation Ti Base / Ceramic Structure</b>
<b>52</b>	<b>Seating and Aftercare</b> Sterilization Intraoral preparation Seating the hybrid abutment and separate crown Seating the hybrid abutment crown Care notes – Implant care
<b>60</b>	<b>General Information</b> Frequently asked questions Material selection table Crystallization and firing parameters Clinical cases


# IPS e.max System

IPS e.max is an innovative all-ceramic system which covers the entire all-ceramic indication range – from thin veneers to 14-unit bridges.

IPS e.max delivers high-strength and highly esthetic materials for the Press and the CAD/CAM technologies. The system consists of innovative lithium disilicate glass-ceramics for smaller restorations and high-strength zirconium oxide for large-span bridges.

The requirements and aims of every case differ. IPS e.max meets these requirements, because due to the system components, you obtain exactly what you need.

- In the field of the **Press technology**, the highly esthetic **IPS e.max Press lithium disilicate glass-ceramic is available** and with **IPS e.max ZirPress a fluorapatite glass-ceramic ingot** for the quick and efficient press-on technique on zirconium oxide.
- For the **CAD/CAM technology**, depending on the case requirements, the innovative **lithium disilicate block IPS e.max CAD** is used or the high-strength **zirconium oxide IPS e.max ZirCAD**.
- The **nano-fluorapatite layering ceramic IPS e.max Ceram**, which is used to characterize and/or veneer all IPS e.max components – glass or oxide ceramics –, completes the IPS e.max System.







# IPS e.max CAD Solutions

IPS e.max CAD stands for individuality. Depending on the indication, users may select from three approaches: This ensures maximum flexibility in the digital work process.

## IPS e.max CAD Monolithic Solutions

Efficient fabrication of full-contour restorations with high strength ( $\geq 360$  MPa) ranging from thin veneers to three-unit bridges.



## IPS e.max CAD Veneering Solutions

Digitally fabricated high-strength veneering structures for zirconium oxide frameworks ( $ZrO_2$ ) – for tooth- and implant-retained crowns and long-span bridges (CAD-on).



**NEW**

## IPS e.max CAD Abutment Solutions

Individual CAD/CAM-fabricated hybrid restorations for implants – for single-tooth restorations in the anterior and posterior region.





# IPS e.max<sup>®</sup> CAD Abutment Solutions

## Product Information

### Description

IPS e.max<sup>®</sup> CAD Abutment Solutions are CAD/CAM-fabricated, implant-supported hybrid restorations for single teeth. These hybrid restorations are individually fabricated of lithium disilicate glass-ceramics (LS<sub>2</sub>) and cemented onto a titanium base (Ti base).

Two different approaches are available:

- IPS e.max CAD hybrid abutment and separate IPS e.max CAD crown
- IPS e.max CAD hybrid abutment crown


Both solutions show outstanding function, efficiency and esthetics. The durable bond to the Ti base is achieved by means of the self-curing Multilink<sup>®</sup> Hybrid Abutment luting composite.

#### Hybrid abutment

The hybrid abutment is an individually milled LS<sub>2</sub> abutment which is luted to the Ti base. The shape, emergence profile and esthetic properties of this abutment can be ideally adjusted to the clinical situation.

Given the lifelike appearance of LS<sub>2</sub> glass-ceramics, the esthetic possibilities are virtually limitless, particularly in the anterior region. Due to the individual characterization, a lifelike appearance is achieved near the root and the transition area to the crown. With the preparation margin of the crown located on the gingival level, the geometry of the hybrid abutments allows for an easy integration of the restoration. Excess cementation material is therefore easily removed.

The milled and crystallized LS<sub>2</sub> ceramic structure is extraorally luted to a Ti base with Multilink Hybrid Abutment, then screwed into place in the oral cavity and finally provided with a permanent IPS e.max CAD crown. Given the convenient fabrication of the hybrid abutment, the process is time-saving and flexible.




#### Hybrid abutment crown

Hybrid abutment crowns are characterized by combining abutment and monolithic crown in one piece. This is an efficient two-in-one solution made of lithium disilicate (LS<sub>2</sub>), which is directly luted to a Ti base.

LS<sub>2</sub> glass-ceramics provide for strength, durability and efficiency, particularly in the posterior region. Moreover, the material offers well-known esthetic properties, allowing restorations to be simply characterized with IPS e.max Ceram stains.

The monolithically milled hybrid abutment crown is reliably and extraorally luted to the Ti base by means of Multilink Hybrid Abutment. Then, the restoration is screwed onto the implant – in one piece. Subsequently, the screw access channel is sealed with a composite material (e.g. Tetric EvoCeram<sup>®</sup>). If required, the screw can be accessed at any time, which affords the dental team clinical flexibility.



IPS e.max CAD hybrid abutment crowns are a new, economically attractive alternative to conventional implant-supported restorations, particularly for the posterior region, where strength, durability and convenient clinical handling matter.

#### Ideally coordinated – Multilink<sup>®</sup> Hybrid Abutment

The self-curing luting composite Multilink Hybrid Abutment in conjunction with Monobond<sup>®</sup> Plus is used for the permanent cementation of ceramic structures made of lithium disilicate glass-ceramic (LS<sub>2</sub>) or zirconium oxide (ZrO<sub>2</sub>) onto bases (e.g. abutment or adhesive base) of titanium/titanium alloy.


This allows

- reliable adhesion due to high bonding values
- optimum esthetics due to two available degrees of opacity
- easy handling due to the convenient Automix syringe

## Material

### IPS e.max CAD

IPS e.max CAD is a lithium disilicate glass-ceramic block for the CAD/CAM technology. It is fabricated using an innovative process which provides an impressive homogeneity of the material. The block can be processed very easily in a CAD/CAM unit in this crystalline intermediate stage. The typical and striking colour of IPS e.max CAD ranges from whitish to blue and bluish-grey. This shade is a result of the composition and the microstructure of the glass-ceramic. The strength of the material in this processable intermediate phase is  $\geq 130$  MPa. After the IPS e.max CAD blocks are milled, the restoration is crystallized in an Ivoclar Vivadent ceramic furnace (e.g. Programat® P500). Unlike with some other CAD/CAM ceramics, the easy-to-conduct crystallization process neither causes any major shrinkage, nor are any complicated infiltration processes required. The crystallization process leads to a change in the microstructure in the IPS e.max CAD material, during which lithium disilicate crystals grow. The densification of 0.2% is accounted for in the CAD software and taken into account upon milling. The final physical properties, such as the strength of  $\geq 360$  MPa and the corresponding optical properties, are achieved through the transformation of the microstructure.



CTE (100–500°C) [ $10^{-6}/K$ ]	$10.5 \pm 0.5$
Flexural strength (biaxial) [MPa]	$\geq 360$ according to ISO 6872
Fracture toughness [ $MPa m^{0.5}$ ]	$\geq 2.0$ according to ISO 6872
Chem. solubility [ $\mu g/cm^2$ ]	$\leq 50$ according to ISO 6872

### Ti base

Ti bases are used for the fabrication of IPS e.max CAD Abutment Solutions. The suitable Ti bases are selected in accordance with the CAD/CAM system used. Please observe the instructions for use and processing of the manufacturer of the Ti bases used.

More information about the authorized CAD/CAM systems is available on the Internet from [www.ivoclarvivadent.com](http://www.ivoclarvivadent.com).



## Uses

### Indications

- Hybrid abutments for anterior and posterior single-tooth restorations
- Hybrid abutment crowns for anterior and posterior restorations

### Contraindications

- Failure to observe the requirements stipulated by the implant manufacturer for using the selected implant type (diameter and length of the implant must be approved for the respective position in the jaw by the implant manufacturer)
- Failure to observe the permissible maximum and minimum ceramic wall thicknesses
- Parafunctions (e.g. bruxism)
- Use of a luting composite other than Multilink® Hybrid Abutment to lute IPS e.max CAD to the Ti base
- Intraoral adhesion of the ceramic structures to the Ti base
- Temporary cementation of the crown on the hybrid abutment
- All uses not stated as indications are contraindicated.

### Important processing restrictions

- Do not mill the blocks with non-compatible CAD/CAM systems.
- If hybrid abutment crowns are fabricated, the opening of the screw channel must not be located in the area of contact points. If this is not possible, a hybrid abutment with a separate crown should be preferred.
- Combination with materials other than IPS e.max Ceram or IPS e.max CAD Crystall./ materials
- Crystallization in a non-recommended ceramic furnace
- Crystallization in a non-calibrated ceramic furnace
- Crystallization in a high-temperature furnace (e.g. Programat® S1)
- Crystallization with deviating firing parameters
- Failure to observe the manufacturer's instructions regarding the processing of the Ti base.

### Side effects

If the patient is known to be allergic to any of the components, IPS e.max CAD and the other materials necessary for the fabrication should not be used.

### Composition

- **IPS e.max CAD blocks**  
Components: SiO<sub>2</sub>, Li<sub>2</sub>O, K<sub>2</sub>O, MgO, Al<sub>2</sub>O<sub>3</sub>, P<sub>2</sub>O<sub>5</sub> and other oxides
- **IPS e.max CAD Crystall./Glaze, Shades and Stains**  
Components: Oxide, glycols
- **IPS e.max CAD Crystall./Glaze Liquid**  
Components: Butandiol
- **IPS e.max CAD Crystall./Add-On**  
Components: Oxides
- **IPS e.max CAD Crystall./Add-On Liquid**  
Components: Water, propylene glycol, butandiol and chloride
- **IPS Object Fix Putty/Flow**  
Components: Oxides, water, thickening agent
- **IPS Natural Die Material**  
Components: Polymethacrylate, paraffin oil, SiO<sub>2</sub> and copolymer
- **IPS Natural Die Material Separator**  
Components: Wax dissolved in hexane
- **Virtual Extra Light Body Fast Set**  
Components: Vinyl polysiloxane, methyl hydrogen polysiloxane, organic platinum complex, silicate and food colouring
- **Multilink Hybrid Abutment**  
Components: Dimethacrylate, HEMA, as well as fillers (barium glass, ytterbium trifluoride, spheroid mixed oxide and titanium dioxide).
- **Monobond Plus**  
Components: Alcohol solution of silane methacrylate, phosphoric acid methacrylate and sulphide methacrylate
- **IPS Ceramic Etching Gel**  
Components: Hydrofluoric acid (approx. 5%)

### Warning

- Do not inhale ceramic dust during finishing. Use exhaust air discharge and mouth protection.
- IPS Ceramic Etching Gel contains hydrofluoric acid. Contact with skin, eyes and clothing must be prevented at all costs, since the material is extremely toxic and corrosive. The etching gel is intended for extraoral use only and must not be applied intraorally (inside the mouth).

## Scientific data

Further scientific data (i.e. strength, wear, biocompatibility) are contained in the **Scientific Documentation IPS e.max CAD**.

The **IPS e.max Scientific Report** contains all studies (in vitro, in vivo) on IPS e.max CAD and the IPS e.max system.

For further information about all-ceramics in general, please refer to the **Ivoclar Vivadent Report No. 16**.




## CAD/CAM partners

IPS e.max CAD has to be processed with an authorized CAD/CAM system. For questions regarding the different systems, please contact the respective cooperation partners.

More information is available on the Internet from [www.ivoclarvivadent.com](http://www.ivoclarvivadent.com).

# IPS e.max CAD Abutment Solutions

Fabrication of IPS e.max CAD hybrid abutment and hybrid abutment crown



## Ivoclar Vivadent Products

Cervitec® Plus, Cervitec® Liquid, Telio® System

OptraGate®, Virtual®

IPS e.max® CAD

Virtual® Extra Light Body Fast Set

IPS e.max® CAD Crystall./ ...  
IPS e.max® Ceram  
Programat® furnaces

IPS Ceramic Etching Gel, Monobond® Plus,  
Multilink® Hybrid Abutment, Liquid Strip

SpeedCEM®, Bluephase®, Tetric EvoCeram®

OptraFine,  
Implant Care


## Shade – tooth shade and abutment shade

Optimum integration in the oral cavity of the patient is the prerequisite for a true-to-nature all-ceramic restoration. To achieve this, the following guidelines and notes must be observed.

With IPS e.max CAD Abutment Solutions, you can imitate not only the clinical crown of a natural tooth, but also a part of the root. This allows you to achieve highly esthetic implant-supported restorations which retain their lifelike appearance also in the case of gingiva recession.


For IPS e.max CAD **hybrid abutment and the separate crown**, the desired tooth shade results from

- the shade of the IPS e.max CAD hybrid abutment (IPS e.max CAD MO ceramic structure, Multilink Hybrid Abutment)
- the shade of the luting material for intraoral cementation of the crown on the IPS e.max CAD hybrid abutment (e.g. SpeedCEM)
- the shade of the IPS e.max CAD LT crown.



For the IPS e.max CAD **hybrid abutment crown**, the desired tooth shade results from

- the shade of the IPS e.max CAD LT ceramic structure
- the shade of Multilink Hybrid Abutment.



## Preparation for the CAD/CAM process

### Scanning

For the fabrication of IPS e.max CAD Abutment Solutions and depending on the CAD/CAM system used, the clinical situation is digitalized either by a direct intraoral scan or an indirect model scan. For notes regarding the scan, please observe the manufacturer's instructions for use of the CAD/CAM system.

### Selecting a Ti base

The required Ti base is selected depending on the inserted implant and the CAD/CAM system used.




## Layer thicknesses of the ceramic components

Observing the geometry requirements of the IPS e.max CAD ceramic structures is the key to success for a durable restoration. The more attention is given to the design, the better the final results and the clinical success will turn out to be.


The following basic guidelines have to be observed:

### Hybrid abutment



- The **wall thickness  $W_{HA}$**  must be at least 0.5 mm.
- The hybrid abutment should be designed in a similar way as a prepared natural tooth:
  - Circular epi-/supragingival shoulder with rounded inner edges or a chamfer.
  - In order for the crown to be cemented to the hybrid abutment using a conventional/self-adhesive cementation protocol, retentive surfaces and a sufficient "preparation height" must be observed.
  - Create an emergence profile with a right angle at the transition to the crown (see picture).
- The **crown width  $B_{Crown}$**  is limited to 6.0 mm from the axial height of contour to the screw channel of the hybrid abutment.
- The notes of the implant manufacturer must be observed regarding the maximum height of the hybrid abutment and separate crown.

### Hybrid abutment crown



- The wall thickness of **hybrid abutment crown  $W_{HAK}$**  must be larger than 1.5 mm for the entire circumference.
- The opening of the screw channel must not be located in the area of contact points. If this is not possible, a hybrid abutment with a separate crown would be preferred.
- The width of the **hybrid abutment crown  $B_{HAK}$**  is limited to 6.0 mm from the axial height of contour to the screw channel.
- The notes of the implant manufacturer must be observed regarding the maximum height of the hybrid abutment crown.



## Block selection

An IPS e.max CAD MO or LT block is selected depending on the indication.

When using a Ti base from Sirona, the dimensions of the interface to the Ti Base (S or L) have to be observed.

<p>IPS e.max CAD hybrid abutment</p>	 <p>IPS e.max CAD <b>MO</b> (Medium Opacity)</p>	<p>Please refer to the table on <b>page 62</b> for the <b>selection of the block shade</b> for the desired tooth shade.</p>
<p>IPS e.max CAD crown (on IPS e.max CAD hybrid abutment)</p>	 <p>IPS e.max CAD <b>LT</b> (Low Translucency)</p>	
<p>IPS e.max CAD hybrid abutment crown</p>	 <p>IPS e.max CAD <b>LT</b> (Low Translucency)</p>	

## CAD/CAM processing

As densification of about 0.2% occurs in IPS e.max CAD during crystallization, this factor has been taken into account in the respective software of the tested CAD/CAM system. Consequently, the milled IPS e.max CAD restorations demonstrate a high accuracy of fit after crystallization. The fabrication steps are described in the directions for use and user manuals of the different CAD/CAM systems. The instructions of the manufacturers must be followed.

## Finishing

It is of critical importance to use the correct grinding instruments for finishing and adjusting the IPS e.max CAD ceramic structure. If unsuitable grinding instruments are used, chipping of the edges and local overheating may occur (please observe the Ivoclar Vivadent Flow Chart "Recommended grinding tools for PS e.max glass-ceramics").

### Basic notes regarding the finishing of IPS e.max CAD

- Carry out adjustments by grinding of IPS e.max CAD restorations while they are still in their pre-crystallized (blue) state, if possible.
- Only use suitable grinding instruments, low speed and light pressure to prevent delamination and chipping at the margins in particular. **Overheating of the glass-ceramic must be avoided.**
- During finishing, make sure that the minimum layer thicknesses are observed.
- Cut the ceramic structure from the block using a diamond separating disc. Slightly scratch the attachment area at the incisal side of the abutment and separate the attachment point from the basal.

### Checking the fit of the ceramic structure on the Ti base

- Carefully place the ceramic structure on the Ti base and check the fit. Observe the position of the rotation lock.



The incisal side of the attachment point is scratched with a diamond separating disc.



The attachment is cut from the basal using a diamond separating disc.



The ceramic structure is carefully placed on the Ti base to check the fit.



Optimum fit of the ceramic structure on the Ti base

## Finishing

### Important!


- Do not finish the shoulder of the ceramic structure to prevent negatively affecting the Ti base.
- Finish the emergence profile if required taking the fit to the gingiva and the minimum thickness (0.5mm) into account.

### Finishing the outer surface of the ceramic structure (hybrid abutment)

- Smooth out the attachment point to the block with fine diamond grinding instruments taking the shape of the emergence profile and the crown margin into account.
- Do not perform any individual shape adjustments, as this will negatively affect the fit of the crown on the hybrid abutment. **Note regarding the crown:** If there are any inaccuracies of fit to the hybrid abutment, adjust the crown by grinding.

### Finishing the outer surface of the ceramic structure (hybrid abutment crown)

- Smooth out the attachment point to the block with fine diamond grinding instruments taking the shape of the emergence profile and the proximal contacts into account.
  - Surface-grind the entire occlusal surface with a fine diamond to smooth out the surface structure created by the CAD/CAM procedure.
  - Check the proximal and occlusal contacts.
  - Design surface textures.
- Clean the ceramic structures with ultrasound in a water bath or blast with the steam jet before further processing.
  - Make sure to thoroughly remove any residue of the milling additive of the CAD/CAM milling unit. Residue of the milling additive remaining on the surface may result in bonding problems and discolouration.
  - Do **not** blast ceramic structures with  $Al_2O_3$  or glass polishing beads!



The shoulder of the ceramic structure must not be finished to prevent negatively affecting the Ti base.




Care should be taken when finishing the emergence profile to prevent affecting the fit to the gingiva.



The attachment point to the block is smoothed out taking the shape of the emergence profile and the crown margin into account.



Individual shape adjustments must not be performed, as this negatively affects the fit of the crown on the hybrid abutment.



The attachment point to the block is smoothed out taking the shape of the emergence profile and the crown margin into account.



The surface of the ceramic structure is ground with a fine diamond to smooth out the surface structure created by the CAD/CAM procedure.

### Tip

Place the crown on the ceramic structure to finish the crown margins. In this way, a smooth transition between the crown and hybrid abutment can be achieved.



next working step ...



**Optional:** Clinical try-in page 17

Completing the IPS e.max CAD ceramic structure page 22

# IPS e.max CAD Abutment Solutions

## Optional: Clinical try-in

A clinical try-in to check the accuracy of fit can be conducted prior to further processing. Clinical try-in may also take place at a later stage, i.e. with the crystallized, tooth-coloured IPS e.max CAD ceramic structure.

### Provisional securing of the ceramic structure on the Ti base

To facilitate the intraoral handling, the components are temporarily attached to one another with silicone material, e.g. Virtual Extra Light Body Fast Set.

The following procedure should be observed in the temporary attachment of the components:

- Clean the untreated Ti base and the ceramic structure with steam and subsequently dry with blown air.
- Place the ceramic structure on the Ti base (which is screwed on the model analog) and mark the relative position of the components with a waterproof pen. This step makes it easier to attain the correct position when the parts are temporarily assembled.
- Seal the screw channel with a foam pellet.
- Insert the Virtual cartridge in the dispenser and remove the protective cap.
- Screw on the mixing tip and attach the Oral Tip to the mixing tip.
- Apply Virtual Extra Light Body Fast Set to the Ti base and directly into the ceramic structure.
- Introduce the Ti base into the ceramic structure. The alignment of the two component must be checked (rotation lock/markings).
- Hold the components firmly in the correct position for 2:30 minutes until Virtual Extra Light Body Fast Set has set.
- Carefully remove any excess that has been displaced with a suitable instrument, e.g. a scalpel.



Cleaned, untreated ceramic structures



The ceramic structure is placed on the Ti base and the relative position is marked.




The screw channel is sealed with a foam pellet.




The Virtual cartridge is inserted in the dispenser. The mixing tip is screwed on and the Oral Tip is attached.



Virtual Extra Light Body Fast Set is applied to the Ti base ...



... and directly on the ceramic structure.




The Ti base is introduced into the ceramic structure. In doing so, the alignment of the two components is checked (rotation lock/markings). The components are firmly held in place for approx. 2:30 minutes until the Virtual Extra Light Body Fast Set has set.



Excess Virtual Extra Light Body Fast Set is removed from the screw channel with an instrument, e.g. a scalpel.



Excess Virtual Extra Light Body Fast Set material is removed from the screw channel with an instrument.



Prepared hybrid abutment or hybrid abutment crown




## Clinical try-in

### Hybrid abutment and dedicated crown

**Important note:** Any intraoral inspection of the occlusion/articulation and necessary grinding adjustments may only be carried out if the components have been attached to one another with Virtual Extra Light Body Fast Set. Virtual has a cushioning effect during the try-in procedure, in particular, if any grinding adjustments have to be made. Therefore, it prevents chipping in the transition area between the hybrid abutment and the crown.

The following procedure should be observed during the clinical try-in:

- The prepared hybrid abutment (provisionally secured in place) and the clean corresponding crown are laid out.
- Remove the provisional restoration.
- Screw the hybrid abutment in manually with the dedicated screw.
- Check the geometry of the hybrid abutment (e.g. fit gingival anaemia) in relation to the gingival margin.
- If desired, the screw channel on the hybrid abutment can be sealed with a foam pellet.
- **Tip:** Isolate the inner aspect of the crown with glycerine gel, e.g. Try-in paste, Liquid Strip
- Place the crown on the hybrid abutment intraorally to check and adjust the proximal contacts, if necessary.
- **Note: No occlusal functional inspection must be performed at this stage.**
- For the functional inspection, the crown has to be secured on the hybrid abutment with Virtual Extra Light Body Fast Set. Try-in paste must not be used for this purpose, as this material is not sufficiently resistant to compressive force.
- Insert the Virtual cartridge in the dispenser and remove the protective cap.
- Screw on the mixing tip and attach the Oral Tip to the mixing tip.
- Apply Virtual Extra Light Body Fast Set to the inner aspect of the crown.
- Press the crown onto the hybrid abutment using the fingers until it has reached the final position. Hold the crown in the final position until the Virtual material has set.
- Remove excess Virtual material.
- Check the occlusion/articulation and make adjustments with suitable grinding instruments, if necessary (see separate IPS e.max recommended grinding instruments for ceramics – use in the dental practice).
- Carefully remove the crown from the hybrid abutment and the hybrid abutment from the implant (including the Ti base).
- Rinse the implant site e.g. with Cervitec Liquid (antibacterial mouth wash with chlorhexidine) to clean and disinfect it.
- Place the temporary restoration.




The hybrid abutment is manually screwed in place with the dedicated screw. The geometry of the hybrid abutment (e.g. fit, gingival anaemia) is checked in relation to the gingival margin.



If desired, the screw channel of the hybrid abutment can be sealed with a foam pellet.



Tip: The inner aspect of the crown can be isolated with glycerine gel.



The crown is placed on the hybrid abutment intraorally to check and if necessary adjust the proximal contacts. **Note: No occlusal functional inspection must be performed at this stage.**



Virtual Extra Light Body Fast Set is applied to the inner aspect of the crown.



The crown is pressed onto the hybrid abutment using the fingers until the final position is reached. The crown is held in the final position until the Virtual material has set.



Excess Virtual material is removed.



The occlusion/articulation is checked and adjustments are made with suitable grinding instruments, if necessary.



The crown is carefully lifted from the hybrid abutment and the Virtual Extra Light Body Fast Set material is removed.



The hybrid abutment is unscrewed.



### Hybrid abutment crown

The following procedure should be observed during the clinical try-in:

- The prepared and cleaned hybrid abutment crown (provisionally secured with in place with Virtual Extra Light Body Fast Set) is laid out.
- Remove the provisional restoration.
- Place the hybrid abutment crown on the implant intraorally in order to check and adjust the proximal contacts, if necessary.
- **Note: No occlusal functional inspection must be performed at this stage.**
- Screw the hybrid abutment crown in manually with the dedicated screw.
- Check the geometry of the hybrid abutment crown (e.g. fit, gingival anaemia) in relation to the gingiva.
- Check the occlusion/articulation and make adjustments with suitable grinding instruments, if necessary (see separate IPS e.max recommended grinding instruments for ceramics – use in the dental practice).
- Carefully remove the hybrid abutment crown.
- Rinse the implant site, e.g. with Cervitec Liquid (antibacterial mouth rinse containing chlorhexidine), to clean and disinfect it.
- Place the temporary restoration.




The hybrid abutment crown is placed on the implant intraorally in order to check and if necessary adjust the proximal contacts. **Note: No occlusal functional inspection must be performed at this stage.**




The hybrid abutment crown is screwed in with the dedicated screw.



The geometry of the hybrid abutment crown is checked (e.g. fit, gingival anaemia) in relation to the gingiva.



The occlusion/articulation is checked and if necessary adjustments are made with suitable grinding instruments.




The hybrid abutment crown (including base) is carefully removed.

# IPS e.max CAD Abutment Solutions

## Completing the IPS e.max CAD ceramic structure

Depending on the desired processing technique and materials, the way to complete the ceramic structure is selected. Basically, two ways to complete the ceramic structure can be distinguished.

- **Polishing technique**  
Polishing of the "blue" restoration, followed by **crystallization** without individual characterization and glaze.
- **Staining technique on the "blue" restoration**  
Characterization and glaze with **IPS e.max CAD Crystall./ materials** on the blue restoration, followed by **Combination firing** (Crystallization and Characterization/Glaze firing in one step).
- **Staining technique on the tooth-coloured restoration**  
**Crystallization** without the application of materials. **Characterization/Glaze firing** of the tooth-coloured restorations with either **IPS e.max CAD Crystall./** or **IPS e.max Ceram materials**.



## Polishing technique



Polishing of the "blue" restoration, followed by crystallization without individual characterization and glaze

If **no characterizations** and **no Glaze firing** are desired, it is possible to polish the ceramic structure manually, followed by crystallization. Please note that polishing causes slight abrasion.

**The polishing technique is preferably used for the emergence profile of the hybrid abutment. For the hybrid abutment crown, the application of glaze is recommended.**



On hybrid abutments, only the emergence profile is polished.



On hybrid abutment crowns, the entire outer aspect is polished.

### Polishing

Please observe the following procedure for polishing the pre-crystallized (blue) ceramic structure:


- Clean the ceramic structure with ultrasound in a water bath or a steam cleaner to remove any contaminations and grease residue.
- Screw Ti base onto a model analog for easier handling.
- Secure the ceramic structure on the Ti base. **Note:** Do **not** finish the Ti base.
- **Overheating of the glass-ceramic must be avoided during polishing.** Observe the recommendations of the manufacturer of the grinding tools.
- Pre-polishing with a diamond rubber polisher (e.g. OptraFine F).
- Fine polishing with a high-gloss rubber polisher (e.g. OptraFine P)
- High-gloss polishing with brushes and polishing paste (e.g. OptaFine HP)
- Clean the ceramic structure with ultrasound in a water bath or the steam jet.



Pre-polishing by means of diamond rubber polishers



Fine polishing by means of high-gloss rubber polishers



High-gloss polishing with brushes and polishing paste



Residue is removed with ultrasound in a water bath...



...or with the steam jet.

### Crystallization

The following steps must be observed:

- Clean the ceramic structure to remove any contaminations and grease residue. Any contamination after cleaning must be prevented.
- Slightly overfill the interface of the ceramic structure with IPS Object Fix Putty or Flow.  
**Immediately reseal the IPS Object Fix Putty/Flow syringe after extruding the material.**
- Place the ceramic structure in the centre of the IPS e.max CAD Crystallization Tray.

### Important

- Conduct the **crystallization** on the IPS e.max CAD Crystallization Tray using the stipulated firing parameters.



**Observe the firing parameters for IPS e.max CAD MO and IPS e.max CAD LT. Firing parameters see page 64**

- **Note:**  
 If a restoration made of IPS e.max CAD MO and one made of IPS e.max CAD LT are to be crystallized in the same firing, the firing parameters for IPS e.max CAD MO must be used.

- Remove ceramic structure from the furnace after completion of the firing cycle (wait for the acoustic signal of the furnace).
- Allow the objects to cool to room temperature in a place protected from draft.
- Do not touch the hot objects with metal tongs.
- Remove the ceramic structure from the IPS e.max CAD Crystallization Tray.
- Remove any residue with ultrasound in a water bath and/or with steam.
- Do **not** remove residue with  $Al_2O_3$  or glass polishing beads.
- Place the ceramic structure on the Ti base and check the fit.
- **If adjustments by grinding of the restoration are required, make sure that no overheating of the ceramic occurs.**



The interface of the ceramic structure is slightly overfilled with IPS Object Fix Putty or Flow. Then the ceramic structure is placed in the centre on the IPS e.max CAD Crystallization Tray.



The crystallization tray is removed from the furnace once the crystallization program has been completed and the object is allowed to cool




The ceramic structure is removed from the IPS e.max CAD Crystallization Tray.



Residue is removed with ultrasound in a water bath....



... or with the steam jet.



Residue **must not** be removed with  $Al_2O_3$  or glass polishing beads.



Polished, crystallized ceramic structure



next working step ...



Permanent cementation Ti base / ceramic structure page 46



## Staining technique on the "blue restoration"



Characterization and glaze with IPS e.max CAD Crystall./ materials on the "blue" restoration, followed by Combination firing

The following paragraphs will explain the steps of glazing and characterizing with IPS e.max CAD Crystall./Shades, Stains and Glaze. In this processing technique, Crystallization and Glaze firing are performed in one step. Characterizations are applied using IPS e.max CAD Crystall./Shades and Stains.



If hybrid abutments are fabricated, only the area of the emergence profile is characterized with IPS e.max CAD Crystall./Shades, Stains and Glaze.



If hybrid abutment crowns are fabricated, the entire outer surface may be individually characterized.

### Required materials

- IPS e.max CAD Crystall./Shades are ready-to-use "Dentin" stains in syringes.



- IPS e.max CAD Crystall./Stains are ready-to-use intensive stains in syringes.



- IPS e.max CAD Crystall./Glaze Paste is a ready-to-use glazing paste.
- IPS e.max CAD Crystall./Glaze Liquid is a special liquid for mixing with Shades, Stains and Glaze.



### Note:

The IPS e.max CAD Crystall./Glaze **Spray** is **not** recommended for glazing IPS e.max CAD Abutment Solutions, as it requires very targeted application. The glazing material must neither reach the bonding surface to the Ti base nor the screw channel, as this may compromise the accuracy of fit.

**Preparing for Combination firing (Crystallization and Stain/Glaze firing in one step)**

- Clean the ceramic structure with the steam jet to remove any contaminations and grease residue. Any contamination after cleaning must be prevented.
- Use the **IPS e.max CAD Crystallization Pin XS** for the crystallization of the ceramic structure.
- Fill the interface of the ceramic structure with either IPS Object Fix Putty or Flow auxiliary firing paste. **Immediately reseal the IPS Object Fix Putty/Flow syringe after extruding the material.**
- Press the **IPS e.max CAD Crystallization Pin XS** **only slightly** into the IPS Object Fix Putty/Flow. **Important: Do not press the pin in too deep to make sure that it does not touch the walls. This may lead to cracks in the ceramic structure.**
- Smooth out displaced auxiliary firing paste using a plastic spatula so that the pin is securely in place.
- Prevent contamination of the outer surface / occlusal surface of the ceramic structure.
- Clean off any possible contamination with a brush dampened with water and dry.




The IPS e.max CAD Crystallization Pin XS is used for the crystallization of the ceramic structure.



The interface of the ceramic structure is filled with either IPS Object Fix Putty or Flow auxiliary firing paste.



**Important:** – The **IPS e.max CAD Crystallization Pin XS** should be pressed **only slightly** into the IPS Object Fix Putty/Flow so that it does not touch the walls of the ceramic structure.



**Incorrect:** Pin pressed in too deep. Pin touches the ceramic structure, which may lead to cracks.



Displaced auxiliary firing paste is smoothed out with a plastic spatula from the margin towards the support pin so that the pin is secured in the paste.



Any possible residue adhering to the outer surface/occlusal surface is cleaned off with a brush dampened with water and dried.



**Combination firing (Crystallization and Stain/Glaze firing in one step)**

Please observe the following procedure for the combination firing:

- Extrude the ready-to-use IPS e.max CAD Crystall./Glaze Paste from the syringe and mix.
- If a slight thinning is desired, the ready-to-use glaze may be mixed with a small amount of IPS e.max CAD Crystall./Glaze Liquid.

**Important:**

- The glazing material must neither reach the bonding surface to the Ti base nor the screw channel, as this may compromise the accuracy of fit. Check the interface before firing and carefully remove any contamination.
- On the abutment, do not apply any materials to the bonding surface to the crown, as this might compromise the fit of the crown.

- Apply IPS e.max CAD Crystall./Glaze Paste evenly on the areas to be glazed using a small brush. Avoid to apply too thick a glaze layer. Avoid "pooling", especially on the occlusal surface of the abutment crown.
- Too thin a glaze layer may lead to an unsatisfactory gloss.
- Apply characterizations with IPS e.max CAD Crystall./Shades and/or IPS e.max CAD Crystall./Stains. For that purpose, extrude the Shades and Stains from the syringe and mix. If necessary, slightly thin them using IPS e.max CAD Crystall./Glaze Liquid. However, the consistency should still remain pasty.
- Apply mixed Shades and Stains directly into the unfired glaze layer using a fine brush. More intensive shades are achieved by several staining procedures and repeated firing, not by applying thicker layers.
- To imitate the incisal area and translucency of the hybrid abutment crown in the incisal and occlusal third, IPS e.max CAD Crystall./Shades Incisal may be used. The cusps and fissures can be individualized using Stains.

**Optional:**

For minor shape adjustments (e.g. proximal or occlusal contact points), IPS e.max CAD Crystall./Add-On is available. The detailed procedure is described on page 33.



After glazing and staining, the Combination firing is conducted in a compatible ceramic furnace (e.g. Programat® CS or Programat P500). When placing the objects into the furnace and setting the firing parameters, observe the following points:

- Place the restoration in the centre of the IPS e.max CAD Crystallization Tray.
- A maximum of 6 units can be positioned on the firing tray and crystallized in the Combination firing with IPS e.max CAD Crystall./Glaze Paste.

**Important**

- Conduct the **Combination firing** on the IPS e.max CAD Crystallization Tray using the stipulated firing parameters.



**Observe the firing parameters for IPS e.max CAD MO and IPS e.max CAD LT. Firing parameters see page 64**

- **Note:**  
If a restoration made of IPS e.max CAD MO and one made of IPS e.max CAD LT are to be crystallized in the same firing, the firing parameters for IPS e.max CAD MO must be used!

- Remove restoration from the furnace after completion of the firing cycle (wait for the acoustic signal of the furnace).
- Allow the objects to cool to room temperature in a place protected from draft.
- Do not touch the hot objects with metal tongs.




IPS e.max CAD Crystall./Glaze Paste is extruded from the syringe and mixed. If required, the paste can be thinned with IPS e.max CAD Crystall./Glaze Liquid.



IPS e.max CAD Crystall./Glaze Paste is applied evenly on the emergence profile of the hybrid abutment or the outer surface of the hybrid abutment crown.



**Important:** The glazing material must reach neither the bonding surface to the Ti base nor the screw channel, as this may compromise the accuracy of fit.



**Important:** The materials must not be applied to the bonding surface to the crown, as this might compromise the fit of the crown.



Individual characterizations of the emergence profile are applied using IPS e.max CAD Crystall./Shades.




Enhancing the chroma on the buccal surface with IPS e.max CAD Crystall./Shades



IPS e.max CAD Crystall./Shade Incisal is applied to imitate the incisal area.



**Optional:** For minor shape adjustments (e.g. proximal contact points), IPS e.max CAD Crystall./Add-On is available.



The ceramic structure is placed in the centre of the IPS e.max CAD Crystallization Tray. The **Combination firing** is conducted using the stipulated **firing parameters**. The firing parameters for IPS e.max CAD MO and IPS e.max CAD LT must be observed.



The ceramic structure is removed from the furnace after completion of the firing cycle (wait for the acoustic signal of the furnace).

## Optional

### Corrective firing

If characterizations or adjustments are required after crystallization, a corrective firing using IPS e.max CAD Crystall./Shades and Stains and Glaze can be conducted. Conduct the corrective firing also on the IPS e.max CAD Crystallization Tray.

For minor shape adjustments (e.g. proximal or occlusal contact points), IPS e.max CAD Crystall./Add-On is available. The detailed procedure is described on page 33.



Once the IPS e.max CAD ceramic structure has cooled to room temperature, proceed with the following steps:

- Remove the ceramic structure from the IPS e.max CAD Crystallization Pin XS.
- Remove any residue with ultrasound in a water bath and/or with the steam jet.
- Do **not** remove residue with Al<sub>2</sub>O<sub>3</sub> or glass polishing beads.
- Place the ceramic structure on the Ti base and check the fit.
- **If adjustments by grinding are required, make sure that no overheating of the ceramic occurs.**




The ceramic structure is removed from the IPS e.max CAD Crystallization Pin XS.



Residue is removed with ultrasound in a water bath. ...



... or with the steam jet.



Residue **must not** be removed with Al<sub>2</sub>O<sub>3</sub> or glass polishing beads.



Glazed and characterized ceramic structures (hybrid abutment crown and hybrid abutment)

## Optional

### Shape adjustments with IPS e.max CAD Crystall./Add-On

For minor shape adjustments (e.g. proximal contact points), IPS e.max CAD Crystall./Add-On is available. The adjustments may be made with both the Combination firing or a separate Corrective firing.

#### Processing

- Mix IPS e.max CAD Crystall./Add-On with IPS e.max CAD Crystall./Add-On Liquid to an easy-to-contour consistency.
- Ensure even mixing of the add-on material and the liquid in order to achieve an optimum firing result.
- Apply the mixed add-on material directly on the unfired Glaze Paste and/or Shades and Stains in the areas to be adjusted and fire.
- Conduct the Combination firing if Add-On is applied on the "blue" partially crystallized restoration.
- Conduct the Corrective firing if Add-On is applied on an already crystallized restoration.



Mixing IPS e.max CAD Crystall./Add-On with IPS e.max CAD Crystall./Add-On Liquid to an easy-to-contour consistency



Application of the mixed Add-On on the blue restoration before crystallization or on the crystallized restoration



Firing parameters see page 64




next working step ...



Permanent cementation Ti base / ceramic structure page 46

## Staining technique on the "tooth-coloured" restoration



Crystallization without application of any materials; separate Characterization/Glaze firing with either IPS e.max CAD Crystall./ or IPS e.max Ceram materials.

### Crystallization

The following steps must be observed:

- Use the **IPS e.max CAD Crystallization Pin XS** for the crystallization of the ceramic structure.
- Fill the interface of the ceramic structure with either IPS Object Fix Putty or Flow auxiliary firing paste. **Immediately reseal the IPS Object Fix Putty/Flow syringe after extruding the material.**
- Slightly press the **IPS e.max CAD Crystallization Pin XS** into the IPS Object Fix Putty/Flow. **Important: Do not press the pin in too deep to make sure that it does not touch the walls. This may lead to cracks in the ceramic structure.**
- Smooth out displaced auxiliary firing paste using a plastic spatula so that the pin is securely in place.
- Prevent contamination of the outer restoration surface. Clean off contamination with a brush dampened with water and dry.
- Place the ceramic structure in the centre of the IPS e.max CAD Crystallization Tray.

### Important

- Conduct the **crystallization** on the IPS e.max CAD Crystallization Tray using the stipulated firing parameters.



Observe the firing parameters for IPS e.max CAD MO and IPS e.max CAD LT. Firing parameters see page 64

### Note:

If a restoration made of IPS e.max CAD MO and one made of IPS e.max CAD LT are to be crystallized in the same firing, the firing parameters for IPS e.max CAD MO must be used!



The IPS e.max CAD Crystallization Pin XS should be used for the crystallization of the ceramic structure.




The interface of the ceramic structure is filled with either IPS Object Fix Putty or Flow auxiliary firing paste.





**Important:** – The IPS e.max CAD Crystallization Pin XS is only slightly pressed into the IPS Object Fix Putty/Flow so that it **does not touch the walls** of the ceramic structure.




**Incorrect:** Pin pressed in too deep. Pin touches the ceramic structure. This may lead to cracks in the ceramic structure.



Displaced auxiliary firing paste is smoothed out with a plastic spatula from the margin towards the support pin so that the pin is secured in the paste.



Any possible residue adhering to the outer surface is cleaned off with a brush dampened with water and dried.



Conduct the crystallization using the stipulated firing parameters. The **firing parameters** for IPS e.max CAD MO and IPS e.max CAD LT must be observed.



Crystallized ceramic structures

- Remove ceramic structures from the furnace after completion of the firing cycle (wait for the acoustic signal of the furnace).
- Allow the objects to cool to room temperature in a place protected from draft.
- Do not touch the hot objects with metal tongs.
- Remove the ceramic structure from the IPS e.max CAD Crystallization Pin.
- Remove any residue with ultrasound in a water bath and/or with the steam jet.
- Do **not** remove residue with Al<sub>2</sub>O<sub>3</sub> or glass polishing beads.
- Place the ceramic structure on the Ti base and check the fit.
- **If adjustments by grinding of the restoration are required, make sure that no overheating of the ceramic occurs.**

next working step, either...

? Stain / Glaze firing with **IPS e.max CAD Crystall.**; page 36

? Stain / Glaze firing with **IPS e.max CAD Ceram**; page 40

### Characterization/Glaze firing with IPS e.max CAD Crystall./...

The following paragraphs will explain the steps of characterizing and glazing with IPS e.max CAD Crystall./Shades, Stains and Glaze.




If hybrid abutments are fabricated, only the area of the emergence profile is characterized with IPS e.max CAD Crystall./Shades, Stains and Glaze.




If hybrid abutment crowns are fabricated, the entire outer surface may be individually characterized.

### Required materials

- IPS e.max CAD Crystall./Shades are ready-to-use "Dentin" stains in syringes.



- IPS e.max CAD Crystall./Stains are ready-to-use intensive stains in syringes.



- IPS e.max CAD Crystall./Glaze Paste is a ready-to-use glazing paste.
- IPS e.max CAD Crystall./Glaze Liquid is a special liquid for mixing with Shades, Stains and Glaze.



#### Note:

The IPS e.max CAD Crystall./Glaze **Spray** is **not** recommended for glazing IPS e.max CAD Abutment Solutions, as it requires very targeted application. The glazing material must neither reach the bonding surface to the Ti base nor the screw channel, as this may compromise the accuracy of fit.



Please observe the following procedure for the Characterization/Glaze firing:

- Extrude the ready-to-use IPS e.max CAD Crystall./Glaze Paste from the syringe and mix.
- If a slight thinning is desired, the ready-to-use glaze may be mixed with a small amount of IPS e.max CAD Crystall./Glaze Liquid.

**Important:**

- The glazing material must neither reach the bonding surface to the Ti base nor the screw channel, as this may compromise the accuracy of fit. Check the interface before firing and carefully remove any contamination.
- On the hybrid abutment, do not apply any materials to the bonding surface to the crown, as this might compromise the fit of the crown.

- Apply IPS e.max CAD Crystall./Glaze Paste evenly on the areas to be glazed using a small brush. Avoid to apply too thick a glaze layer. Avoid "pooling", especially on the occlusal surface of the hybrid abutment crown.
- Too thin a glaze layer may lead to an unsatisfactory gloss.
- Apply characterizations with IPS e.max CAD Crystall./Shades and/or IPS e.max CAD Crystall./Stains. For that purpose, extrude the Shades and Stains from the syringe and mix. If necessary, slightly thin them using IPS e.max CAD Crystall./Glaze Liquid. However, the consistency should still remain pasty.
- Apply mixed Shades and Stains directly into the unfired glaze layer using a fine brush. More intensive shades are achieved by several staining procedures and repeated firing, not by applying thicker layers.
- To imitate the incisal area and translucency of the abutment crown in the incisal and occlusal third, IPS e.max CAD Crystall./Shades Incisal may be used. The cusps and fissures can be individualized using Stains.

After glazing and staining, the Characterization/Glaze firing (Corrective firing) is conducted in a compatible ceramic furnace (e.g. Programat® CS or Programat P500). When placing the objects into the furnace and setting the firing parameters, observe the following points:

- Place the restoration in the centre of the IPS e.max CAD Crystallization Tray.
- A maximum of 6 units can be positioned on the firing tray for the firing with IPS e.max CAD Crystall./Glaze Paste.

**Important**

- Conduct the **Corrective firing** on the IPS e.max CAD Crystallization Tray using the stipulated firing parameters.




Firing parameters see page 64




IPS e.max CAD Crystall./Glaze Paste is extruded from the syringe and mixed. If required, the paste is thinned with IPS e.max CAD Crystall./Glaze Liquid.




IPS e.max CAD Crystall./Glaze Paste is applied evenly on the emergence profile of the hybrid abutment or the outer surface of the hybrid abutment crown.



**Important:** The glazing material must reach neither the bonding surface to the Ti base nor the screw channel or the bonding surface to the crown, as this may compromise the accuracy of fit.



Characterizing the emergence profile with Shades



Enhancing the chroma of the buccal surface



IPS e.max CAD Crystall./Shade Incisal is applied to imitate the incisal area.



The Corrective firing is conducted on the IPS e.max CAD Crystallization Tray using the stipulated firing parameters.

## Optional

### Corrective firing

- If adjustments are required, another Corrective firing using IPS e.max CAD Crystall./Shades and Stains and Glaze can be conducted. Conduct the Corrective firing also on the IPS e.max CAD Crystallization Tray.
- For minor shape adjustments (e.g. proximal contact points), IPS e.max CAD Crystall./Add-On is available. The adjustments may be made with both Crystallization/Glaze and Corrective firing.
- The detailed procedure is described on page 33.



Once the IPS e.max CAD ceramic structure has cooled to room temperature, proceed with the following steps:


- Remove the ceramic structure from the IPS e.max CAD Crystallization Pin XS.
- Remove any residue with ultrasound in a water bath and/or with the steam jet.
- Do **not** remove residue with  $Al_2O_3$  or glass polishing beads.
- Place the ceramic structure on the Ti base and check the fit.
- If adjustments by grinding are required, make sure that no overheating of the ceramic occurs.
- If the restoration is ground, manually polish the corresponding areas to a high gloss after grinding.



The ceramic structure is removed from the IPS e.max CAD Crystallization Pin XS.



Any residue is removed with ultrasound in a water bath or with the steam jet.



Residue **must not** be removed with  $Al_2O_3$  or glass polishing beads.



Glazed and characterized ceramic structures (hybrid abutment and/or hybrid abutment crown)



next working step ...



Permanent cementation Ti base / ceramic structure page 46

## Characterization/Glaze firing with IPS e.max Ceram

The following paragraphs will explain the steps of characterizing and glazing with **IPS e.max Ceram**.



If hybrid abutments are fabricated, only the area of the emergence profile is characterized with IPS e.max Ceram Shades, Essences and Glaze.



If hybrid abutment crowns are fabricated, the entire outer surface may be individually characterized with IPS e.max Ceram Shades, Essences, and Glaze.

### Required materials

- **IPS e.max Ceram Essences** are intensively shaded stains in powder form.
- **IPS e.max Ceram Shades** are ready-to-use stains in syringes.
- **IPS e.max Ceram Glaze and Stain Liquid (allround, longlife)** to mix the materials in powder form (Essences, Glaze), as well as to thin paste materials (Shades, Glaze).



**IPS e.max CAD Crystall./Shades, Stains, Glaze and IPS e.max Ceram Shades, Essence, Glaze must not be mixed with each other!**

The following steps must be observed:

- Clean the finished ceramic structure with the steam jet to remove any contaminations and grease residue. Any contamination after cleaning must be prevented.
- For better wetting of the stains, a small quantity of IPS e.max Ceram Glaze and Stain Liquid may be slightly rubbed into the area that needs to be characterized.
- Mix the pastes or powders with the IPS e.max Ceram Glaze and Stain Liquid allround or longlife to the desired consistency.
- More intensive shades are achieved by several staining procedures and repeated firing, not by applying thicker layers.
- To imitate the incisal area and translucency of the hybrid abutment crown in the incisal and occlusal third, IPS e.max Ceram Shade Incisal may be used. The cusps and fissures can be individualized using Essences.
- If hybrid abutments are fabricated, only the area of the emergence profile is characterized with IPS e.max Ceram Shades and Essences.
- Secure the ceramic structure on the firing pin of the honey-comb tray with a little IPS Object Fix Putty or Flow for firing.


### Important:

The characterization must neither reach the bonding surface to the Ti base nor the screw channel, as this may compromise the accuracy of fit. Check the interface before firing and carefully remove any contamination. On the hybrid abutment, do not apply any materials to the bonding surface to the crown, as this may compromise the fit of the crown.




**Conduct the Characterization/Glaze firing for IPS e.max Ceram on a honey-comb firing tray using the stipulated firing parameters. Firing parameters see page 64**






IPS e.max Ceram Shade Incisal is applied to imitate the incisal area.



Enhancing the chroma of the buccal surface



Individual characterization of the emergence profile with IPS e.max Ceram Essences



The Stain and Characterization firing is conducted on a honey-comb firing tray.

- Remove restoration from the furnace after completion of the firing cycle (wait for the acoustic signal of the furnace).
- Allow the objects to cool to room temperature in a place protected from draft.
- Do not touch the hot objects with metal tongs.

Additional Characterization firings can be conducted with the same firing parameters.

#### Glaze firing

Glaze firing is conducted with powder or paste glaze. On abutments, only the emergence profile is glazed. On hybrid abutment crowns, glaze is applied to the entire outer surface.

#### Required materials

- **IPS e.max Ceram Glaze Paste, Glaze Powder** are glazing materials in paste and powder forms.
- **IPS e.max Ceram Glaze and Stain Liquid (allround, longlife)** to mix the materials in powder form (Essences, Glaze), as well as to thin paste materials (Shades, Glaze)



**IPS e.max CAD Crystall./Shades, Stains, Glaze and IPS e.max Ceram Shades, Essence, Glaze must not be mixed with each other!**



The following procedure is recommended:

- For easier handling, the ceramic structure can be positioned on the Ti base for glazing. For that purpose, secure Ti base on a model analog.
- Mix the glazing material (IPS e.max Ceram Glaze Paste or Powder) with the IPS e.max Ceram Glaze and Stain Liquid allround or longlife to the desired consistency.
- Apply an even layer of glazing material covering all areas that are to be glazed.
- If required, the fluorescence may be increased by applying a fluorescing glazing material (paste or powder).

**Important:**

The glazing material **must neither reach the bonding surface** to the Ti base nor the screw channel, as this may compromise the accuracy of fit. Check the interface before firing and carefully remove any contamination. On the abutment, do not apply any glaze to the bonding surface to the crown, as this might compromise the fit of the crown.




The **IPS e.max CAD Crystall./Glaze Spray** is not recommended for glazing IPS e.max CAD Abutment Solutions, as it requires very targeted application. The glazing material must neither reach the bonding surface to the Ti base nor the screw channel, as this may compromise the accuracy of fit.




Conduct the Characterization/Glaze firing for IPS e.max Ceram on a honey-comb firing tray using the stipulated firing parameters. Firing parameters see page 64



An even layer of glaze material is applied to the emergence profile of the hybrid abutment. Care has to be taken that no glaze material enters the screw channel.



The glazing material is applied evenly on the outer surface of the hybrid abutment crown. Care has to be taken that no glaze material enters the screw channel.



Care has to be taken that no glaze material is present on the interface of the hybrid abutment and hybrid abutment crown prior to the firing cycle. The glaze material is carefully removed, if necessary.



The Characterization/Glaze firing is conducted on a honey-comb firing tray with the corresponding parameters.


- Remove restoration from the furnace after completion of the firing cycle (wait for the acoustic signal of the furnace).
- Allow the objects to cool to room temperature in a place protected from draft.
- Do not touch the hot objects with metal tongs.

## Optional

### Shape adjustments of IPS e.max Ceram Add-On

Use IPS e.max Ceram Add-On Dentin and/or Incisal for shape adjustments after Glaze firing. Please observe the following procedure for processing:

- Mix IPS e.max Ceram Add-On Dentin or Incisal with IPS e.max Ceram Build-Up Liquid soft or allround and apply on the corresponding areas.
- Fire with the stipulated parameters for the "Add-On after Glaze firing". Observe long-term cooling!
- If necessary, polish the adjusted areas to a high gloss after firing.



Firing parameters see page 64



next working step ...



Permanent cementation Ti base / ceramic structure page 46

# IPS e.max<sup>®</sup> CAD Abutment Solutions

## Crown on IPS e.max CAD hybrid abutment


The crown on the IPS e.max Hybrid Abutment can be completed using either the staining technique or the cut-back technique. To characterize and glaze, either the PS e.max CAD Crystall./ materials or the IPS e.max Ceram materials are used.

**Basically, the procedure for completing a crown is the same as that for a crown on a prepared tooth. For detailed information about the procedure, please refer to the IPS e.max CAD Instructions for Use.**

**Example: IPS e.max CAD crown – Cut-back technique – IPS e.max Ceram**



Partially reduced IPS e.max CAD restorations fitted on the model. Always observe minimum thicknesses!



For crystallization, the partially reduced IPS e.max CAD restorations are placed directly on the IPS e.max CAD Crystallization Tray using IPS Object Fix Putty or Flow.



The wash firing is conducted using e.g. IPS e.max Ceram Glaze, Shades and Essences.





Completion of the anatomical shape of the reduced areas using IPS e.max Ceram Incisal and Opal materials



Finishing with diamond burs and design of a true-to-nature shape and surface structure Finally, glaze firing is conducted using IPS e.max Ceram Glaze.



IPS e.max CAD crown after glaze firing (partially reduced and veneered with IPS e.max Ceram) on a IPS e.max CAD hybrid abutment

# IPS e.max CAD Abutment Solutions

## Permanent cementation of base / ceramic structure

Careful preparation of the bonding surface is a prerequisite for the successful adhesive cementation of the base and the ceramic structure. The following paragraphs outline the required procedure. The procedure is the same for hybrid abutments and hybrid abutment crowns.

### Required materials

- IPS Ceramic Etching Gel
- Monobond® Plus
- Multilink® Hybrid Abutment
- Glycerigel (z.B. Liquid Strip)



	IPS e.max CAD ceramic structure (LS <sub>2</sub> )	Base
<b>Blasting</b>	-	According to the instructions of the manufacturer
<b>Etching</b>	Bonding area to the base with IPS® Ceramic Etching Gel for 20 s	-
<b>Conditioning</b>	The bonding area with Monobond® Plus for 60 s	
<b>Adhesive cementation</b>	Multilink® Hybrid Abutment	
<b>Covering the cementation joint</b>	Glycerine gel, e.g. Liquid Strip	
<b>Curing</b>	7 minutes auto-polymerization	
<b>Polishing the cementation joint</b>	Conventional polishers for ceramic/composite resin	

### Preparation of the Ti base

The following procedure should be observed in the preparation of the Ti base for the cementation with the ceramic structure:

- The Ti base should be prepared according to the instructions of the manufacturer.
- Clean the Ti base with an ultrasonic bath or with a steam cleaner and then dry it with blown air.
- Screw the Ti base on the model analog.
- Place the ceramic structure on the Ti base and mark the relative position of the components with a waterproof pen. This facilitates locating the correct position when the parts are assembled at a later stage.
- The emergence profile of the base must not be blasted or modified in any way.
- **If the manufacturer recommends that the bonding surface of the Ti base be blasted, the following procedure should be observed:**
  - Apply silicone (Virtual Extra Light Body Fast Set) to protect the emergence profile and the screw channel .
  - Carefully blast the bonding area according to the instructions of the manufacturer.
  - Remove silicone.
  - Clean the Ti base with ultrasound in a water bath or with the steam jet.
  - After cleaning, the bonding surface must not be contaminated under any circumstances, as this would impair the bond.
- Apply Monobond Plus on the clean bonding surface and allow it to react for 60 s. After the reaction time, disperse any residue with air that is free of water and oil.
- Seal the screw channel with a foam pellet or wax. The bonding surface must not be contaminated in the process.



The Ti base is screwed on the model analog. The relative position to the ceramic structure is marked with a waterproof pen.



E.g. silicone (Virtual Extra Light Body Fast Set) is applied in order to protect the emergence profile and the screw channel.



The bonding surface can be carefully blasted **according to the instructions of the manufacturer.**



Removal of the silicone and subsequently cleaning with ultrasound in a water bath or with the steam jet.



Monobond Plus is applied to the clean bonding surface and allowed to react for 60 s. After the reaction time, any remaining residue is dried with blown air that is free of water and oil.




The screw channel is sealed with a foam pellet or wax. The bonding surface must not be contaminated in the process.

### Preparing the ceramic structure

The following procedure must be observed in the preparation of the ceramic structure for cementation on the Ti base:

- Do **not** blast the ceramic structure in preparation for the cementation.
- Clean the ceramic structure in an ultrasonic bath or with a steam cleaner and subsequently blow dry.
- After cleaning, the bonding surface must not be contaminated under any circumstances, as this would impair the bond.
- Wax can be applied to protect the outer surfaces or the glazed areas.
- Etch the bonding surface with 5% hydrofluoric acid gel (IPS Ceramic Etching Gel) for 20 s.
- Subsequently, thoroughly rinse the bonding surface under running water and dry with air that is free of water and oil.
- Apply Monobond Plus on the clean bonding surface and allow it to react for 60 s. After the reaction time, dry any remaining residue with blown air that is free of water and oil.



The ceramic structure **must not** be blasted.



Etching with IPS Ceramic Etching Gel for 20 seconds. Subsequently, the restoration is rinsed with water and blown dry.



Monobond Plus is allowed to react for 60 s, and excess is blown dry.

### Cementation with Multilink® Hybrid Abutment

The following instructions must be observed in the cementation procedure:

- The cleaned and conditioned components (ceramic structure, Ti base) are laid out ready for cementation.
- **The subsequent cementation procedure must be carried out quickly and without interruption. The working time of Multilink Hybrid Abutment is approximately 2 min at 23 °C (± 1 °C) or 73 °F (± 1.8 °F).**
- As a general rule, a new mixing tip is attached to the Multilink Hybrid Abutment syringe prior to each use.
- Apply a thin layer of Multilink Hybrid Abutment directly from the mixing syringe to the bonding surface of the Ti base and **the bonding surface of the ceramic structure.**
- The mixing tip is left on the Multilink Hybrid Abutment syringe until the next use. The remaining cement polymerizes in the tip and functions as a seal.
- Place the ceramic structure on the Ti base in such a way that the position markings are aligned.
- Press the parts lightly and evenly together and check the correct relative position of the components (transition Ti base/ ceramic structure).
- Subsequently, press the parts tightly together for 5 s.
- Carefully remove excess in the screw channel, e.g. with a Microbrush or brush, using rotary movements.

#### Important:

- **Important: Excess must not be removed before curing has started, i.e. 2-3 minutes after mixing. For the purpose, a suitable dental lab instrument (e.g. Le Cron) is used. The components are held in place with light pressure in the process.**
- Glycerine gel is applied (e.g. Liquid Strip) on the cementation joint to prevent the formation of an inhibition layer. The glycerine gel must be applied cautiously to avoid blending it with or displacing the composite. Make sure to leave the gel on the cementation joint until polymerization is complete.
- Next, the luting composite is completely auto-polymerized within 7 min.
- **Important: Do not move the components until Multilink Hybrid Abutment has completely cured. They can be held immobile with e.g. diamond-coated tweezers.**
- After the completion of auto-polymerization, rinse off the glycerine gel with water.
- **Make sure to cautiously polish the cementation joint with rubber polishers at a low speed (< 5,000 rpm) to avoid overheating.**
- Remove any cement residue left in the screw channel with suitable rotating instruments.
- Clean the restoration with ultrasound in a water bath or with the steam jet.




The cleaned and conditioned components are laid out ready for cementation.



A new mixing tips is attached to the Multilink Hybrid Abutment prior to each use. The Multilink Hybrid Abutment mixing syringe is attached.



A thin layer of Multilink Hybrid Abutment is directly applied from the mixing tip to the bonding surface of the Ti base.




A thin layer of Multilink Hybrid Abutment is directly applied from the mixing tip on the bonding surface of the ceramic structure.




The ceramic structure is placed on the Ti base in such a way that the position markings are aligned. The components are joined using even and light pressure and the relative position of the components is checked (transition base/ceramic structure).



Subsequently, the components are tightly pressed together for 5 s.




Excess in the screw channel is carefully removed, e.g. with a Microbrush or brush, using rotary movements.



**Important:** Excess must not be removed before curing has started, i.e. 2-3 minutes after mixing. The components are held in place with light pressure in the process.



Glycerine gel (e.g. Liquid Strip) is applied on the cementation joint to prevent the formation of an inhibition layer.



The luting composite auto-polymerizes within 7 min. **Important:** The components must not be moved until auto-polymerization is completed. The components must be immobilized during this time.



After the completion of auto-polymerization, the glycerine gel is rinsed off with water.



The cementation joint is cautiously polished with rubber polishers at low speed (< 5,000 rpm), to avoid overheating.



Any remaining cement residue in the screw channel is removed with suitable rotating instruments. The Ti base must not be damaged.



Completed hybrid abutment and hybrid abutment crown after cementation



# IPS e.max CAD Abutment Solutions

## Seating and Aftercare

### Sterilization

The hybrid abutments or hybrid abutment crowns must be sterilized prior to insertion. Furthermore, the locally applicable legal regulations and the hygiene standards applicable for a dental practice must be observed.

Steam sterilization can be performed with the 3 x fractionated pre-vacuum with the following parameters: Sterilization time 3 min; steam temperature 132 °C/270 °F; resulting in a half-cycle exposure time of 1.5 min. The abutment is for immediate use. No storage after sterilization!



The responsibility for the sterility of the hybrid abutment or hybrid abutment crown lies with the user. It must be ensured that only suitable devices, materials and product-specifically validated methods are used to perform sterilization. The equipment and devices must be properly maintained and serviced at regular intervals. The fabricator (dental technician) of the IPS e.max CAD Abutment Solution must inform the dentist of the need to sterilize the abutment before inserting it in the patient's mouth!

### Intraoral preparation

Please observe the following procedure to prepare for the permanent cementation of the implant-supported restoration:

- Remove the temporary restoration.
- Clean the implant site.
- Check the periimplant tissue (emergence profile).



## Seating the hybrid abutment and crown

### Preparing/conditioning the hybrid abutment with dedicated crown

Conditioning of the ceramic surface, i.e. the bonding surface, in preparation for cementation is critical for generating a sound bond between the cementation material and the all-ceramic material.

The following procedure must be observed in the preparation of the ceramic structure for cementation on the Ti base:

- Do **not** blast IPS e.max CAD hybrid abutment or IPS e.max CAD crown with  $\text{Al}_2\text{O}_3$  or glass polishing beads.
- Ideally, the clinical try-in is conducted before etching to prevent contamination of the bonding surface.
- Thoroughly clean the hybrid abutment and crown with water and subsequently blow dry.
- Etch the bonding surface with 5% hydrofluoric acid gel (IPS Ceramic Etching Gel) for 20 s. The etching gel must not come into contact with the emergence profile or the outer side of the crown. **Important: No intraoral application of the IPS Ceramic Etching Gel.**
- Subsequently, thoroughly rinse the bonding surface under running water and dry it with air that is free of water and oil.
- If an adhesive or self-adhesive cementation protocol is used, apply Monobond Plus to the clean bonding surface and allow it to react for 60 s. After this reaction time, disperse any residue with air that is free of water and oil.



The IPS e.max CAD ceramic structures **must not** be blasted in preparation for cementation.



The bonding surfaces are etched with IPS Ceramic Etching Gel and subsequently cleaned.



Monobond Plus is applied to the bonding surfaces, and allowed to react for 60 s. Excess is dispersed with air.

## Seating the hybrid abutment and dedicated crown

### Note:

**Temporary insertion of the IPS e.max CAD crown on the IPS e.max CAD hybrid abutment is contraindicated!**

For the permanent seating of the hybrid abutment and the crown, please observe the following working steps. Please also observe the Instructions for Use of the selected luting material.

**SpeedCEM® is recommended for the seating of IPS e.max CAD crowns on IPS e.max hybrid abutments.**



- Do not use phenolic mouth washes, as such products negatively influence the bond between the ceramic and the composite.
- Insert the hybrid abutment intraorally into the implant.
- Manually screw in the matching implant screw.
- Tighten the implant screw with a torque wrench (observe the instructions of the manufacturer).
- Insert a cotton or foam pellet into the screw channel.
- Seal the screw channel with a temporary composite (e.g. Telio® CS Inlay). This serves to ensure access to the screw at a later stage.
- Check the bonding area for contamination/moisture and clean or dry with an air syringe, if necessary.
- Apply the luting material, e.g. **SpeedCEM**, into the conditioned crown.
- Place the crown onto the hybrid abutment and secure in place in the final position.
- Conduct the pre-polymerization using the four-quarter technique.
- Remove excess luting material.
- Cover the cementation joint with glycerine gel (e.g. Liquid Strip).
- Polymerize with an LED curing light (e.g. Bluephase®).
- Rinse off the glycerine gel with water.
- Check the occlusion and articulation and make adjustments, if necessary. If adjustments are made to the restoration by grinding, these areas must subsequently be polished to a high gloss, e.g. using OptraFine.
- Polish restoration margins and the cementation joint with silicone polishers (e.g. Astropol®, OptraFine).
- Apply Cervitec Plus in the area of the gingival margin.



The hybrid abutment is inserted into the implant intraorally.




The matching implant screw is screwed in manually.



The implant screw is tightened with a torque wrench (the instructions of the manufacturer must be observed).



The screw channel is sealed, for instance with a cotton or foam pellet and a temporary composite material.



The luting material, e.g. SpeedCEM, is applied into the conditioned crown.



The crown is placed on the hybrid abutment and secured in place.



Pre-polymerization using the four-quarter technique



Excess luting material is removed.



The restoration margin is covered with glycerine gel (e.g. Liquid Strip).



The luting material is cured with an LED curing light (e.g. Bluephase).



The glycerine gel is rinsed off with water.



The occlusion and articulation is checked and adjustments are made, if necessary.



The restoration margins and the cementation joint are polished (e.g. OptraPol, OptraFine).



Completed IPS e.max CAD hybrid abutment and crown

## Seating the hybrid abutment crown

### Preparing/Conditioning the Hybrid Abutment Crown

Please observe the following notes to prepare for the intraoral sealing of the screw channel:

- As a general rule, do **not** blast IPS e.max CAD hybrid abutment crowns with  $\text{Al}_2\text{O}_3$  or glass polishing beads.
- Thoroughly clean the the hybrid abutment crown with water and blow dry.
- Etch the screw channel from the occlusal side with 5% hydrofluoric acid gel (IPS Ceramic Etching Gel) for 20 seconds. Make sure that no etching gel comes into contact with the occlusal surface. **Important: Do not use the IPS Ceramic Etching Gel intraorally.**
- Thoroughly rinse off the etching gel with water and dry with oil- and water-free air.
- Apply Monobond Plus to the etched and cleaned surface in the screw channel, allow to react for 60 seconds and then disperse excess with oil- and water-free air.



IPS e.max CAD ceramic structures **must not** be blasted.



The screw channel is etched with IPS Ceramic Etching Gel for 20 s and subsequently cleaned.



Monobond Plus is applied, allowed to react for 60 s and excess is dispersed.

### Seating the hybrid abutment crown

For the permanent seating of the hybrid abutment crown, please observe the following working steps:

- Do not use phenolic mouth washes, as such products negatively influence the bond between the ceramic and the composite.
- Insert the hybrid abutment crown intraorally into the implant.
- Manually screw in the matching implant screw.
- Tighten the implant screw with a torque wrench (observe the instructions of the manufacturer).
- Check the screw channel for contamination/moisture and clean with Total Etch (phosphoric acid gel), if necessary.
- Insert a cotton or foam pellet into the screw channel.
- Apply the bonding agent, followed by polymerization.
- Seal the screw channel with a composite material (e.g. Tetric EvoCeram) in the appropriate shade.
- Polymerize with an LED curing light (e.g. Bluephase).
- Check the occlusion/articulation after polymerization and correct possible rough spots with suitable fine-grit diamonds.
- Polish to a high gloss with silicone polishers (e.g. OptraPo/OptraFine).





The hybrid abutment crown is inserted into the implant intraorally.



The matching implant screw is screwed in manually.



The implant screw is tightened with a torque wrench (the instructions of the manufacturer must be observed).



The screw channel is sealed with a composite material (e.g. Tetric EvoCeram) in the appropriate shade.



Polymerization with an LED curing light (e.g. Bluephase)



After polymerization, the occlusion/articulation is checked and possible rough spots are adjusted with suitable finishers or fine diamonds.



High-gloss polishing is performed using silicone polishers (e.g. Astropol P, Astropol HP or Astrobrush).



Completed IPS e.max CAD hybrid abutment crown

## Care notes – Implant Care

Implant Care comprises a range of coordinated products for the professional care of patients during the various phases of implant treatment and lifelong aftercare. Products for professional tooth cleaning and bacterial control contribute to the long-term quality assurance of implant-supported restorations. Structural elements peri-implant tissue, natural teeth, dental restorations, gingiva and the mucosa are treated in an optimum way with regard to function and esthetics.



# IPS e.max<sup>®</sup> CAD Abutment Solutions

## General Information

### Frequently Asked Questions

#### In addition to the desired tooth shade, why should the root shade also be defined/determined upon shade determination?

IPS e.max CAD Abutment Solutions allow you to fabricate restorations with a lifelike appearance both in the visible area and the area below the gingiva (root). By defining the root shade, a highly esthetic outcome can be achieved especially in the case of receding gingiva.

#### Is it possible to fabricate an abutment or an abutment crown with IPS e.max CAD (LS<sub>2</sub>) without using a Ti base?

**No!** For this indication, IPS e.max CAD needs the support provided by the Ti base. In addition, the Ti base allows an optimum (industrially fabricated) fit to the implant to be achieved.

#### Which Ti bases can be used for the fabrication of IPS e.max CAD Abutment Solutions?

Only Ti bases of authorized CAD/CAM systems may be used. More information about the CAD/CAM cooperation systems is available on the Internet from [www.ivoclarivadent.com](http://www.ivoclarivadent.com).

#### Is it permissible to modify the selected Ti base?

The Ti base must not be adjusted by grinding, as this would compromise the fit of the IPS e.max CAD ceramic structure. The instructions of the manufacturer regarding the preparation for permanent cementation must be observed.

#### Is a hybrid abutment crown indicated in the anterior region?

This indication depends on the position and inclination of the implant. If the opening of the screw channel is located on the oral surface, a hybrid abutment crown may be fabricated in the anterior region.

#### May a hybrid abutment crown be cut-back and subsequently supplemented with IPS e.max Ceram layering materials.?

No. For implant-supported restorations, it is recommended to fabricate monolithic restorations (without veneer). In this way, chipping of the layering ceramic is prevented.

#### Do IPS e.max CAD ceramic structures have to be glazed in all cases?

No. High gloss can also be achieved by a corresponding polishing procedure. The polishing technique (before crystallization) is preferably used for the emergence profile of the hybrid abutment. For the hybrid abutment crown, the application of glaze is recommended.

#### Is it possible to use an IPS e.max CAD hybrid abutment as an abutment for a bridge restoration?

No. Only single-tooth restorations may be fabricated.

#### Can different CAM units be used for milling the IPS e.max CAD ceramic structure (abutment) and the dedicated IPS e.max CAD crown?

If different CAM units are used, inaccuracies of fit may occur in unfavourable cases. Therefore, both IPS e.max CAD objects (abutment, crown) should be ideally milled in the same CAM unit.

#### Can a clinical try-in be conducted with the IPS e.max CAD Abutment Solutions? How are the ceramic structures prepared for this?

Yes. Clinical try-in may be performed either before or after crystallization of the IPS e.max CAD ceramic structures. The Ti base and IPS e.max CAD ceramic structure are temporarily joined in the laboratory by means of a silicone material, e.g. Virtual Extra Light Body Fast Set. This facilitates the intraoral handling during clinical try-in with the patient.

#### What must be observed for the clinical try-in of a crown on a hybrid abutment?

To check the occlusion/articulation and to make possible adjustments, the crown must be temporarily secured on the hybrid abutment with a silicone material, e.g. Virtual Extra Light Body Fast Set. The silicone material acts as a buffer and prevents chipping in the marginal area of the crown. Try-in pastes or Vaseline must not be used for functional checks.



**Can a glaze spray be used for glazing IPS e.max CAD ceramic structures (e.g. IPS e.max CAD Crystall./Glaze Spray)?**

*We do not recommend using the Glaze Spray for the indications hybrid abutment or hybrid abutment crown, as there is a risk that the bonding surface to the Ti base or the screw channel are contamination with glaze.*

**What material is used to permanently cement the IPS e.max CAD ceramic structures to the Ti base?**

*Only Multilink Hybrid Abutment is to be used for permanent cementation. This ensures a high-quality bond. Given the high opacity of the luting composite, complete optical masking of the Ti base is achieved and thus an excellent esthetic appearance ensured.*

**How is the Ti base prepared for the permanent cementation with Multilink Hybrid Abutment?**

*Provided it has been approved by the manufacturer, carefully blast the bonding area with  $Al_2O_3$  at low pressure until an even mat surface has been achieved. After cleaning, the area is conditioned with Monobond Plus.*

**How is the screw channel of a hybrid abutment crown sealed after seating?**

*The screw channel is extraorally conditioned (etching, silanating). After the restoration has been intraorally screwed down on the implant, the screw channel is sealed with a restorative composite.*

# Material Selection Table

## IPS e.max CAD hybrid abutment and IPS e.max CAD dedicated crown

The material is selected on the basis of the desired tooth shade (Bleach BL or A-D). Depending on the geometry of the hybrid abutment and the crown, shade adjustment by means of characterization with IPS e.max CAD Crystall/Shades, Stains, or IPS e.max Ceram Shades and Essences may be necessary to achieve the desired shade.

The block recommendations for the hybrid abutment have been selected in such a way that the desired tooth shade is achieved in combination with the crown. In the "cervical area", it may be necessary to characterize the hybrid abutment according to the clinical situation.

Material combination to achieve the tooth shade	Bleach BL and A-D Shade Guide																			
	BL1	BL2	BL3	BL4	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Extraoral cementation IPS e.max CAD abutment / Ti base	Ti base																			
IPS e.max CAD ceramic structure	MultiLink Hybrid Abutment HO 0*																			
Intraoral cementation Crown on hybrid abutment	Adhesive, self-adhesive or conventional cementation, e.g. SpeedCEM																			
IPS e.max CAD crown	MO 0	MO 1	MO 2	MO 3	MO 1	MO 3	MO 3	MO 1	MO 1	MO 3	MO 3	MO 1	MO 1	MO 1	MO 1	MO 1	MO 1	MO 1	MO 1	MO 3
	LT BL1	LT BL2	LT BL3	LT BL4	LT A1	LT A2	LT A3	LT A3.5	LT A4	LT B1	LT B2	LT B3	LT B4	LT C1	LT C2	LT C3	LT C4	LT D2	LT D3	LT D4

\*The range of products may vary from country to country.

## IPS e.max CAD hybrid abutment crown

The material is selected on the basis of the desired tooth shade (Bleach BL or A-D). Depending on the geometry of the hybrid abutment crown, shade adjustment by means of characterization with IPS e.max CAD Crystall/Shades, Stains, or IPS e.max Ceram Shades and Essences may be necessary to achieve the desired shade. In the "cervical area", it may be necessary to characterize the hybrid abutment crown according to the clinical situation.

Material combination to achieve the tooth shade	Bleach BL and A-D Shade Guide																			
	BL1	BL2	BL3	BL4	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Extraoral cementation IPS e.max CAD abutment crown / Ti base	Ti base																			
IPS e.max CAD ceramic structure	MultiLink Hybrid Abutment HO 0*																			
	LT BL1	LT BL2	LT BL3	LT BL4	LT A1	LT A2	LT A3	LT A3.5	LT A4	LT B1	LT B2	LT B3	LT B4	LT C1	LT C2	LT C3	LT C4	LT D2	LT D3	LT D4

\*The range of products may vary from country to country.

1 IPS e.max CAD LT Blocks are available in 10 shades. To create the desired tooth shade, select the closest block shade in the respective shade group and determine the restoration shade by means of Stains.

## Clinical Cases

### IPS e.max CAD hybrid abutment / IPS e.max CAD-crown

Dr. R. Watzke / F. Perkon, Ivoclar Vivadent, Liechtenstein



Starting situation with shaped emergence profile



IPS e.max CAD ceramic structure (abutment) / IPS e.max CAD crown, milled



IPS e.max CAD hybrid abutment / IPS e.max CAD crown, completed



Screwed in IPS e.max CAD hybrid abutment



IPS e.max CAD crown on IPS e.max CAD hybrid abutment, cemented

### IPS e.max CAD hybrid abutment crown

Dr. L. Enggist / P. Scherrer, Ivoclar Vivadent, Liechtenstein



Starting situation



IPS e.max CAD hybrid abutment crowns (prepared for clinical try-in)



Try-in of the IPS e.max CAD hybrid abutment crowns



Completed IPS e.max CAD hybrid abutment crowns



Seated IPS e.max CAD hybrid abutment crowns



## Crystallization and Firing Parameters

### Crystallization/Combination firing: IPS e.max CAD MO – optional for LT

with or without application of IPS e.max CAD Crystall./ materials



Furnaces	Stand-by temperature B [°C/°F]	Closing time S [min]	Heating rate t <sub>1</sub> [°C/°F/min]	Firing temperature T <sub>1</sub> [°C/°F]	Holding time H <sub>1</sub> [min]	Heating rate t <sub>2</sub> [°C/°F/min]	Firing temperature T <sub>2</sub> [°C/°F]	Holding time H <sub>2</sub> [min]	Vacuum 1 T <sub>1</sub> [°C/°F] T <sub>2</sub> [°C/°F]	Vacuum 2 Z <sub>1</sub> [°C/°F] Z <sub>2</sub> [°C/°F]	Long-term cooling L [°C/°F]	Cooling rate t <sub>3</sub> [°C/°F/min]
Programat P300 P500 P700	403/757	06:00	60/108	770/1418	00:10	30/54	850/1562	10:00	550/770 1022/1418	770/850 1418/1562	700/1292	0
Programat CS Program 7	403/757	06:00	60/108	770/1418	00:10	30/54	850/1562	10:00	550/770 1022/1418	770/850 1418/1562	700/1292	0

### Crystallization/Combination firing: IPS e.max CAD LT – not suitable for MO

with or without application of IPS e.max CAD Crystall./ materials



Furnaces	Stand-by temperature B [°C/°F]	Closing time S [min]	Heating rate t <sub>1</sub> [°C/°F/min]	Firing temperature T <sub>1</sub> [°C/°F]	Holding time H <sub>1</sub> [min]	Heating rate t <sub>2</sub> [°C/°F/min]	Firing temperature T <sub>2</sub> [°C/°F]	Holding time H <sub>2</sub> [min]	Vacuum 1 T <sub>1</sub> [°C/°F] T <sub>2</sub> [°C/°F]	Vacuum 2 Z <sub>1</sub> [°C/°F] Z <sub>2</sub> [°C/°F]	Long-term cooling L [°C/°F]	Cooling rate t <sub>3</sub> [°C/°F/min]
Programat P300 P500 P700	403/757	06:00	90/162	820/1508	00:10	30/54	840/1544	07:00	550/820 1022/1508	820/840 1508/1544	700/1292	0
Programat CS Program 1	403/757	06:00	90/162	820/1508	00:10	30/54	840/1544	07:00	550/820 1022/1508	820/840 1508/1544	700/1292	0

### Corrective firing – Characterization/Glaze firing IPS e.max CAD MO, LT

with IPS e.max CAD Crystall./ materials



Furnaces	Stand-by temperature B [°C/°F]	Closing time S [min]	Heating rate t <sub>1</sub> [°C/°F/min]	Firing temperature T <sub>1</sub> [°C/°F]	Holding time H <sub>1</sub> [min]	Heating rate t <sub>2</sub> [°C/°F/min]	Firing temperature T <sub>2</sub> [°C/°F]	Holding time H <sub>2</sub> [min]	Vacuum 1 T <sub>1</sub> [°C/°F] T <sub>2</sub> [°C/°F]	Vacuum 2 Z <sub>1</sub> [°C/°F] Z <sub>2</sub> [°C/°F]	Long-term cooling L [°C/°F]	Cooling rate t <sub>3</sub> [°C/°F/min]
Programat P300 P500 P700	403/757	06:00	90/162	820/1508	00:10	30/54	840/1544	03:00	550/820 1022/1508	820/840 1508/1544	700/1292	0
Programat CS Program 3	403/757	06:00	90/162	820/1508	00:10	30/54	840/1544	03:00	550/820 1022/1508	820/840 1508/1544	700/1292	0

### Characterization/Glaze firing

with IPS e.max Ceram Shades, Essences, Glaze



Furnaces	Stand-by temperature B [°C/°F]	Closing time S [min]	Heating rate t <sub>1</sub> [°C/°F/min]	Firing temperature T <sub>1</sub> [°C/°F]	Holding time H [min]	Vacuum 1 V <sub>1</sub> [°C/°F]	Vacuum 2 V <sub>2</sub> [°C/°F]	Long-term cooling L [°C/°F]
Programat P300 P500 P700	403/757	06:00	60/108	770/1418	1:00 – 2:00	450/842	769/1416	500/932

**Note:**

If the layer thickness is less than 2 mm on the IPS e.max CAD object, long-term cooling (L) is not required.

### Corrective firing

with IPS e.max Ceram Add-On



Furnaces	Stand-by temperature B [°C/°F]	Closing time S [min]	Heating rate t <sub>1</sub> [°C/°F/min]	Firing temperature T <sub>1</sub> [°C/°F]	Holding time H [min]	Vacuum 1 V <sub>1</sub> [°C/°F]	Vacuum 2 V <sub>2</sub> [°C/°F]	Long-term cooling L [°C/°F]
Programat P300 P500 P700	403/757	06:00	50/90	700/1292	01:00	450/842	699/1290	500/932

## Crystallization and Firing Parameters

The following points should be observed for ceramic furnaces used for the crystallization of IPS e.max CAD:

- Crystallization should be carried out in an Ivoclar Vivadent ceramic furnace (e.g. Programat CS, Programat P300, P500, or P700).
- If you use other, untested ceramic furnaces, please consult Ivoclar Vivadent about their compatibility with IPS e.max CAD.

Basically, the following applies:

- Ceramic furnaces without
  - function for controlled (long-term) cooling
  - vacuum functioncannot be used.
- Before the first crystallization and every six months after that, the ceramic furnace must be calibrated.
- Depending on the mode of operation, more frequent calibrations may be required. Observe the instructions of the manufacturer.

The following aspects should be observed for conducting the **crystallization**:

- Use only IPS Object Fix Putty or Flow as an auxiliary firing paste to place the restoration directly on the IPS e.max CAD Crystallization Tray.
- IPS e.max CAD restorations must not be directly placed on the IPS e.max CAD Crystallization Tray and the Pins, i.e. without auxiliary firing paste, for crystallization.
- Use only the IPS e.max CAD Crystallization Tray and the corresponding Pins, since they store the heat necessary for slow and above all tension-free cooling of the glass-ceramic.
- Always conduct the crystallization under vacuum.
- Remove IPS e.max CAD objects from the furnace after completion of the firing cycle (wait for the acoustic signal of the furnace).
- Allow the objects to cool to room temperature in a place protected from draft.
- Do not touch the hot objects with metal tongs.
- Do not blast or quench the objects.



crystallization and firing parameters



please fold out the page



# Ivoclar Vivadent – worldwide

## **Ivoclar Vivadent AG**

Bendererstrasse 2  
9494 Schaan  
Liechtenstein  
Tel. +423 235 35 35  
Fax +423 235 33 60  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

## **Ivoclar Vivadent Pty. Ltd.**

1 – 5 Overseas Drive  
P.O. Box 367  
Noble Park, Vic. 3174  
Australia  
Tel. +61 3 979 595 99  
Fax +61 3 979 596 45  
[www.ivoclarvivadent.com.au](http://www.ivoclarvivadent.com.au)

## **Ivoclar Vivadent Ltda.**

Alameda Caiapós, 723  
Centro Empresarial Tamboré  
CEP 06460-110 Barueri – SP  
Brazil  
Tel. +55 11 2424 7400  
Fax +55 11 3466 0840  
[www.ivoclarvivadent.com.br](http://www.ivoclarvivadent.com.br)

## **Ivoclar Vivadent Inc.**

1-6600 Dixie Road  
Mississauga, Ontario  
L5T 2Y2  
Canada  
Tel. +1 905 670 8499  
Fax +1 905 670 3102  
[www.ivoclarvivadent.us](http://www.ivoclarvivadent.us)

## **Ivoclar Vivadent Shanghai**

**Trading Co., Ltd.**  
2/F Building 1, 881 Wuding Road,  
Jing An District  
200040 Shanghai  
China  
Tel. +86 21 6032 1657  
Fax +86 21 6176 0968  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

## **Ivoclar Vivadent Marketing Ltd.**

Calle 134 No. 7-B-83, Of. 520  
Bogotá  
Colombia  
Tel. +57 1 627 33 99  
Fax +57 1 633 16 63  
[www.ivoclarvivadent.co](http://www.ivoclarvivadent.co)

## **Ivoclar Vivadent SAS**

B.P. 118  
F-74410 Saint-Jorioz  
France  
Tel. +33 450 88 64 00  
Fax +33 450 68 91 52  
[www.ivoclarvivadent.fr](http://www.ivoclarvivadent.fr)

## **Ivoclar Vivadent GmbH**

Dr. Adolf-Schneider-Str. 2  
D-73479 Ellwangen, Jagst  
Germany  
Tel. +49 (0) 79 61 / 8 89-0  
Fax +49 (0) 79 61 / 63 26  
[www.ivoclarvivadent.de](http://www.ivoclarvivadent.de)

## **Wieland Dental + Technik GmbH & Co. KG**

Schwenninger Strasse 13  
D-75179 Pforzheim  
Germany  
Tel. +49 (0) 72 31 / 37 05-0  
Fax +49 (0) 72 31 / 35 79 59  
[www.wieland-dental.com](http://www.wieland-dental.com)

## **Ivoclar Vivadent Marketing (India) Pvt. Ltd.**

503/504 Raheja Plaza  
15 B Shah Industrial Estate  
Veera Desai Road, Andheri (West)  
Mumbai, 400 053  
India  
Tel. +91 (22) 2673 0302  
Fax +91 (22) 2673 0301  
[www.ivoclarvivadent.in](http://www.ivoclarvivadent.in)

## **Ivoclar Vivadent s.r.l.**

Via Isonzo 67/69  
40033 Casalecchio di Reno (BO)  
Italy  
Tel. +39 051 611 35 55  
Fax +39 051 611 35 65  
[www.ivoclarvivadent.it](http://www.ivoclarvivadent.it)

## **Ivoclar Vivadent K.K.**

1-28-24-4F Hongo  
Bunkyo-ku  
Tokyo 113-0033  
Japan  
Tel. +81 3 6903 3535  
Fax +81 3 5844 3657  
[www.ivoclarvivadent.jp](http://www.ivoclarvivadent.jp)

## **Ivoclar Vivadent Ltd.**

12F W-Tower, 1303-37  
Seocho-dong, Seocho-gu,  
Seoul 137-855  
Republic of Korea  
Tel. +82 (2) 536 0714  
Fax +82 (2) 596 0155  
[www.ivoclarvivadent.co.kr](http://www.ivoclarvivadent.co.kr)

## **Ivoclar Vivadent S.A. de C.V.**

Av. Insurgentes Sur No. 863.  
Piso 14, Col. Napoles  
03810 México, D.F.  
México  
Tel. +52 (55) 50 62 10 00  
Fax +52 (55) 50 62 10 29  
[www.ivoclarvivadent.com.mx](http://www.ivoclarvivadent.com.mx)

## **Ivoclar Vivadent BV**

De Fruittuinen 32  
2132 NZ Hoofddorp  
Netherlands  
Tel. +31 23 529 37 91  
Fax +31 23 555 45 04  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

## **Ivoclar Vivadent Ltd.**

12 Omega St, Rosedale  
PO Box 303011 North Harbour  
Auckland 0751  
New Zealand  
Tel. +64 9 914 99 99  
Fax +64 9 914 99 90  
[www.ivoclarvivadent.co.nz](http://www.ivoclarvivadent.co.nz)

## **Ivoclar Vivadent Polska Sp. z o.o.**

Al. Jana Pawla II 78  
00-175 Warszawa  
Poland  
Tel. +48 22 635 54 96  
Fax +48 22 635 54 69  
[www.ivoclarvivadent.pl](http://www.ivoclarvivadent.pl)

## **Ivoclar Vivadent Marketing Ltd.**

Prospekt Andropova 18 korp. 6/  
office 10-06  
115432 Moscow  
Russia  
Tel. +7 499 418-03-00  
Fax +7 499 418-03-10  
[www.ivoclarvivadent.ru](http://www.ivoclarvivadent.ru)

## **Ivoclar Vivadent Marketing Ltd.**

Qlaya Main St.  
Siricon Building No.14, 2<sup>nd</sup> Floor  
Office No. 204  
P.O. Box 300146  
Riyadh 11372  
Saudi Arabia  
Tel. +966 1 293 83 45  
Fax +966 1 293 83 44  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

## **Ivoclar Vivadent Pte. Ltd.**

171 Chin Swee Road  
#02-01 San Centre  
Singapore 169877  
Tel. +65 6535 6775  
Fax +65 6535 4991  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

## **Ivoclar Vivadent S.L.U.**

c/ Emilio Muñoz Nº 15  
Entrada c/ Albaracin  
E-28037 Madrid  
Spain  
Tel. + 34 91 375 78 20  
Fax + 34 91 375 78 38  
[www.ivoclarvivadent.es](http://www.ivoclarvivadent.es)

## **Ivoclar Vivadent AB**

Dalvägen 14  
S-169 56 Solna  
Sweden  
Tel. +46 (0) 8 514 93 930  
Fax +46 (0) 8 514 93 940  
[www.ivoclarvivadent.se](http://www.ivoclarvivadent.se)

## **Ivoclar Vivadent Liaison Office**

: Tesvikiye Mahallesi  
Sakayik Sokak  
Nisantas' Plaza No:38/2  
Kat:5 Daire:24  
34021 Sisli – Istanbul  
Turkey  
Tel. +90 212 343 08 02  
Fax +90 212 343 08 42  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

## **Ivoclar Vivadent Limited**

Ground Floor Compass Building  
Feldspar Close  
Warrens Business Park  
Enderby  
Leicester LE19 4SE  
United Kingdom  
Tel. +44 116 284 78 80  
Fax +44 116 284 78 81  
[www.ivoclarvivadent.co.uk](http://www.ivoclarvivadent.co.uk)

## **Ivoclar Vivadent, Inc.**

175 Pineview Drive  
Amherst, N.Y. 14228  
USA  
Tel. +1 800 533 6825  
Fax +1 716 691 2285  
[www.ivoclarvivadent.us](http://www.ivoclarvivadent.us)



Manufacturer:  
Ivoclar Vivadent AG, 9494 Schaan, Liechtenstein  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

Date information prepared: 2013-04/Rev. 0

Some products and/or indications may not be regulatory cleared/released in all markets. Please contact the local Ivoclar Vivadent sales office for the current national status.

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

Printed in Liechtenstein  
© Ivoclar Vivadent AG, Schaan / Liechtenstein  
662738/en

Traditional 510(k)  
Straumann® Variobase™ Abutments

Straumann USA, LLC  
January 23, 2014

  
**ivoclar**  
**vivadent**  
technical

A17-67

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 18

---

**Appendix 18 – Lava™ Frame Instructions for Use**







**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 19

---

**Appendix 19 – IPS e.max® Press Abutment Solutions 510(k)  
Summary**



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room –WO66-G609  
Silver Spring, MD 20993-0002

Ivoclar Vivadent AG  
C/O Ms. Donna Marie Hartnett  
Director of Quality Assurance / Regulatory Affairs  
Ivoclar Vivadent, Incorporated  
175 Pineview Drive  
Amherst, New York 14228

OCT 18 2012

Re: K120053

Trade/Device Name: IPS e.max<sup>®</sup> Press – Abutment Solutions  
Regulation Number: 21 CFR 872.3630  
Regulation Name: Endosseous Dental Implant Abutment  
Regulatory Class: II  
Product Code: NHA, EIH  
Dated: September 28, 2012  
Received: October 2, 2012

Dear Ms. Hartnett:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Page 2- Ms. Hartnett

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Anthony D. Watson, B.S., M.S., M.B.A.  
Director  
Division of Anesthesiology, General Hospital,  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

K 12 0053

**Indications for Use**

510(k) Number (if known): K120053

Device Name: IPS e.max® Press – Abutment Solutions

**Indications For Use:**

IPS e.max® Press Abutment Solutions is intended for use in partially or fully endentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations.

IPS e.max Press Abutment Solutions is recommended for the fabrication of:

- Hybrid abutments for single-tooth restorations
- Hybrid abutment crowns for restorations


The following Ti bases are intended to be used with IPS e.max Press Abutment Solutions.

Implant manufacture	Implant system, diameter	Compatible Ti base (abutment), dimensions diameter = D gingiva height = GH height = HTi
Straumann®	Bone Level RC Ø 4.1 mm or Ø 4.8 mm 510K K062129	RC Cementable abutment D 5.0–6.5 mm GH 1.0–3.0 mm HTi 4.0–5.5 mm 510K K072071

Prescription Use X AND/OR  
(Part 21 CFR 801 Subpart D)

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)  
  
 \_\_\_\_\_  
 (Division Sign-Off)  
 Division of Anesthesiology, General Hospital  
 Infection Control, Dental Devices

Page 1 of 1

510(k) Number: K120053

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 20

---

**Appendix 20 – P.004 RC Cementable Abutments 510(k) Summary**

K072071

1072

Section H  
510(k) SUMMARY

AUG 15 2007

**1. Applicant's Name and Address**

Straumann US (on behalf of Institut Straumann AG)  
60 Minuteman Rd.  
Andover, MA 01810  
Telephone Number: 978-747-2513  
Fax Number: 978-747-0023  
Contact Person: Elaine Alan  
Regulatory Affairs Specialist  
Date of Submission: July 27, 2007

**2. Name of the Device**

Trade Name: P.004 RC Cementable Abutments  
Common Name: Abutment, Dental, Endosseous implants  
Classification Name: Abutment, Dental, Endosseous implants  
Regulation Number: 21 CFR 872.3630

**3. Legally Marketed Device to which Equivalence is Claimed (Predicate Device)**

P.004 RC Anatomic Abutments, K062129  
Temporary Copings, K041070  
ITI Protective Healing Caps, K962023

**4. Description of the Device**

The Straumann P.004 Dental Implant System is an integrated system of endosseous dental implants, which are designed to support prosthetic devices for partially or fully edentulous patients. The system consists of a variety of dental implants, abutments and surgical and prosthetic parts and instruments.

Abutments are placed into dental implants to provide support for prosthetic restorations, temporary copings support a temporary restoration out of occlusion and protective caps protect the abutment during the healing phase.

**5. Intended Use of the Device**

Abutments are placed into dental implants to provide support for prosthetic restorations such as crowns, bridges and overdentures. Abutments can be used in single tooth replacements and multiple tooth restorations. Temporary Copings are intended to serve as a base for temporary restorations for up to 28 days. Protective Caps are intended to protect the outer configuration of the abutment and to maintain and

510(k) Submission: RC Cementable Abutments  
July 27, 2007

Straumann US  
Page 39

2 of 2  
K072071

condition the contours of the soft tissue during the healing phase for up to 6 months.

**6. Technological Characteristics**

The modified Anatomic Abutments, Temporary Copings and Protective Caps are substantially equivalent to the currently cleared devices. The intended use is **identical** to the predicate devices. The proposed devices have the same material composition, basic design and fundamental operating principles to the currently cleared devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Institut Straumann AG  
C/O Ms. Elaine Alan  
Regulatory Affairs Specialist  
Straumann USA  
60 Minuteman Road  
Andover, Massachusetts 01810

AUG 15 2007

Re: K072071  
Trade/Device Name: P.004 Cementable Abutments, Temporary Copings and  
Protective Caps  
Regulation Number: 872.3630  
Regulation Name: Endosseous Dental Implant Abutment  
Regulatory Class: II  
Product Code: NHA  
Dated: July 27, 2007  
Received: July 30, 2007

Dear Ms. Alan:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.



Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

1071

K072071

**INDICATIONS FOR USE STATEMENT**

Device Name: P.004 Cementable Abutments

Indications for Use:

Abutments are intended to be placed into dental implants to provide support for prosthetic reconstructions such as crowns, bridges and overdentures.

Device Name: P.004 Protective Caps

Indications for Use:

Protective Caps are intended to protect the outer configuration of the abutment and to maintain and condition the contours of the soft tissue during the healing phase for up to 6 months.

Device Name: P.004 Temporary Copings

Indications for Use:


Temporary Copings are intended to serve as a base for temporary restorations for up to 28 days.

Prescription Use   X   AND/OR Over-The-Counter Use         
(Part 21 CFR 801 Subpart D) (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of 1

  
\_\_\_\_\_  
Division Sign-Off  
Division of Anesthesiology, General Hospital,  
Infection Control, Dental Devices

510(k) Submission: RC Cementable Abutments 510(k) Number:   K072071   Straumann US  
July 27, 2007 Page 5

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 21

---

**Appendix 21 – IPS e.max® CAD Abutment Solutions 510(k)**  
**Summary**

**510(K) SUMMARY (revised)**  
**IPS e.max® CAD Abutment Solutions**

Contact: Donna Marie Hartnett, Director QA/Regulatory Affairs

Company: Ivoclar Vivadent, 175 Pineview Drive, Amherst, NY 14228  
(716) 691-0010

Date Prepared: October 18, 2013

Proprietary Name: IPS e.max® CAD Abutment Solutions

Classification Name: Abutment, Implant, Dental Endosseous (872.3630)

Predicate Devices: Sirona Dental CAD/CAM System (K111421)  
 IPS e.max CAD (K051705)  
 IPS e.max Press Abutment Solutions (K120053 and K124008)

OCT 31 2013

**Device Description:** IPS e.max CAD Abutment Solutions is intended for use in partially or fully edentulous mandibles and maxillae in support of single cement-retained restorations. Titanium bases are a premanufactured prosthetic component directly connected to dedicated endosseous dental implants. The Titanium base is used for adhesion to mesostructures to restore function and esthetics in the oral cavity. IPS e.max CAD Abutment Solutions are lithium disilicate blocks in various sizes. One side of the block is mounted to a mandrel that will be inserted into the spindle's clamping chuck of the grinding machine. The connection geometry to titanium bases is prefabricated, i.e. already include in the shipped block. Connection geometries fit select Titanium Bases marketed by Straumann, Nobel Biocare and Biomet 3i as identified in the Intended Use section. The mesostructure is individually designed and milled using CAD/CAM Technology into the shape of a hybrid abutment or hybrid abutment crown as designed by the trained professional using the Sirona inLab and Cerec SW 4.2 (or higher) software.. The device serves as the esthetic mesostructure which is extraorally cemented onto a Titanium Base. The two piece abutment is mounted onto the implant and fixed with a screw.

**Predicate Device:** The predicate device to which IPS e.max® CAD Abutment Solutions has been compared is Sirona Dental CAD/CAM System(K111421). For this application, IPS e.max® CAD Abutment Solutions has been compared to its predicate with regard to chemical composition, performance data and indications for use. The comparison shows that IPS e.max® CAD Abutment Solutions is substantially equivalent to the predicate device.

510K Summary 5-1

## 510(K) SUMMARY (revised)

### IPS e.max® CAD Abutment Solutions

Intended Use:

IPS e.max CAD Abutment Solutions is intended for use in partially or fully edentulous mandibles and maxillae in support of single cement-retained restorations. The system comprises three parts:

IPS e.max CAD ceramic structure,

Ti base and

CAD/CAM software

The IPS e.max CAD ceramic structure cemented to the Ti base is recommended for two-piece hybrid abutments for single tooth restorations and hybrid abutment crowns, used in conjunction with endosseous dental implants. The compatible Implant systems, Ti bases and CAD/CAM systems are shown below:

Implant systems: Nobel Biocare Replace (K020646), Nobel Biocare Activem (K071370), Straumann Bone Level (K053088, ), Biomet 3i Certain (K014235)

CAD/CAM Systems: Sirona inLab and Cerec SW 4.2 (or higher) software

Titanium bases:

Implant manufacturer	Implant System	Implant Diameter (mm)	TiBase	Sirona Ref.	Interface size
Nobel Biocare	Replace NP	3.5	NBRS 3.5	6282474	L
	Replace RP	4.3	NBRS 4.3	6282482	L
	Replace WP	5.0	NBRS 5.0	6282490	L
	Replace 6.0	6.0	NBRS 6.0	6282508	L
Nobel Biocare	Nobel Active NP	3.5	NB A 4.5	6208188	L
	Nobel Active RP	4.3 / 5.0	NB A 5.0	6208253	L
Straumann	Bone Level NC	3.3	S BL 3.3	6308154	L
	Bone Level RC	4.1 / 4.8	S BL 4.1	6308337	L
Biomet 3i	Certain	3.4	B C 3.4	6308048	S
	Certain	4.1	B C 4.1	6308097	L
	Certain	5.0	B C 5.0	6308121	L

For the titanium base Straumann Bone Level 3.3 L the indication is restricted for replacement of single lateral incisors in the maxilla and lateral and central incisors in the mandible.

Material Composition: The device is composed of proprietary lithium disilicate (Li Si<sub>2</sub>) dental ceramic and is identical in composition to IPS e/max CAD (K051705).

## **510(K) SUMMARY (revised)**

IPS e.max® CAD Abutment Solutions

**Technological Characteristics:** The device design, i.e. delivery form, and intended use of IPS e.max CAD Abutment Solutions and the predicate device are the same. The materials comply with ISO 6872:2008 for Dental Ceramics. The composition of the subject device has been modified from the predicate, however, there are no ingredients in the subject device which pose any new issues of safety and effectiveness.

**Scientific Concept:** The underlying scientific concept is the use of an already introduced technology of a titanium base abutment combined with individually CAD/CAM fabricated ceramic prosthetics made from lithium disilicate, a material proven to be suitable for safe and effective dental restoratives.

**Testing Summary:** According to FDA Guidance “Class II Special Controls Guidance Document: Root-form Endosseous Dental Implants and Endosseous Dental Implant Abutments” may 12, 2004, fatigue testing has been performed for angled abutments.

**Physical properties of IPS e.max CAD:**

- CTE (100°C – 500°C)             $10.5 \pm 0.5 \cdot 10^{-6}/K$
- Flexural strength (Biaxial)       $\geq 360 \text{ MPa}$  (Test Method ISO 6872)
- Fracture toughness                 $\geq 2.0 \text{ MPa m}^{0.5}$  (Test Method ISO 6872)
- Chemical solubility                $\leq 50 \text{ } \mu\text{g}/\text{cm}^2$  (Test Method ISO 6872)
- Crystallization temperature       $840 - 850^\circ\text{C}$

**Conclusion:** IPS e.max CAD Abutment Solutions is substantially equivalent to the predicate device.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-002

October 31, 2013

Ivoclar Vivadent AG  
C/O Donna Hartnett  
Director of Quality Affairs/Regulatory Affairs  
Ivoclar Vivadent, Incorporation  
175 Pineview Dr.  
AMHERST NY, 14228

Re: K132209  
Trade/Device Name: IPS e.max® CAD Abutment Solutions  
Regulation Number: 21 CFR 872.3630  
Regulation Name: Endosseous dental implant abutment  
Regulatory Class: II  
Product Code: NHA  
Dated: August 6, 2013  
Received: August 7, 2013

Dear Ms. Hartnett:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act

or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Mary S. Bunner -S

Kwame Ulmer, MS  
Director  
Division of Anesthesiology, General Hospital,  
Respiratory, Infection Control and  
Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure



## Indications for Use

**510(k) Number (if known):**           K132209          

**Device Name:** IPS e.max® CAD Abutment Solutions

**Indications For Use:**

IPS e.max CAD Abutment Solutions is intended for use in partially or fully edentulous mandibles and maxillae in support of single cement-retained restorations. The system comprises three parts:

- IPS e.max CAD mesostructure,
- Ti base and
- CAD/CAM software.

The IPS e.max CAD mesostructure cemented to the Ti base is recommended for two-piece hybrid abutments for single tooth restorations and hybrid abutment crowns, used in conjunction with endosseous dental implants. The compatible Implant systems, Ti bases and CAD/CAM systems are shown below:

**Implant systems:** Nobel Biocare Replace (K020646), Nobel Biocare Activem (K071370), Straumann Bone Level (K053088, K062129, K060958), Biomet 3i Osseotite (K980549)

**CAD/CAM Systems:** Sirona inLab and Cerec SW 4.2 and above

**Titanium bases:**

Implant manufacturer	Implant System	Implant Diameter (mm)	TiBase	Sirona Ref.	Interface size
Nobel Biocare	Replace NP	3.5	NBRS 3.5	6282474	L
	Replace RP	4.3	NBRS 4.3	6282482	L
	Replace WP	5.0	NBRS 5.0	6282490	L
	Replace 6.0	6.0	NBRS 6.0	6282508	L
Nobel Biocare	Nobel Active NP	3.5	NB A 4.5	6208188	L
	Nobel Active RP	4.3 / 5.0	NB A 5.0	6208253	L
Straumann	Bone Level NC	3.3	S BL 3.3	6308154	L
	Bone Level RC	4.1 / 4.8	S BL 4.1	6308337	L
Biomet 3i	Certain	3.4	B C 3.4	6308048	S
	Certain	4.1	B C 4.1	6308097	L
	Certain	5.0	B C 5.0	6308121	L

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of 2

For the titanium base Straumann Bone Level 3.3 L the indication is restricted for replacement of single lateral incisors in the maxilla and lateral and central incisors in the mandible.

Prescription Use X AND/OR  
(Part 21 CFR 801 Subpart D)

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF  
NEEDED)

Andrew L Steen S  
2013.10.30 13:55:24 -04'00'

Concurrence of CDRH, Office of Device Evaluation (ODE)  
Page 2 of 2

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 22

---

**Appendix 22 – Sirona Dental CAD/CAM System 510(k) Summary**

K111421

FEB 17 2012

**510(k) Summary**  
**for**  
**Sirona Dental Systems**  
**Sirona Dental CAD/CAM System**

**1 Sponsor**

Sirona Dental Systems GmbH

Fabrikstrasse 31

D-64625 Bensheim

Germany

Contact Person: Fritz Kolle

Telephone: 49 6251 16 32 94

Date Prepared: May 06, 2011

**2 Device Name**

Proprietary Name: Sirona Dental CAD/CAM-System

Common/Usual Name: Abutment, implant, dental, endosseous

Classification Names: Endosseous dental implant abutment

**3 Predicate Devices**

Replace® NP, K091756, Brånemark®, K091756, Tissue level NN, K081005, OsseoSpeed™, K081666, Frialit® / Xive®, K032158, Osseotite K072642, Tapered Screw-Vent®, K060880, Nobel Active NP, K102436, Bone Level NC, K062129, Certain®, K073345.

Sirona Dental CAD/CAM System (K100152)

---

Sirona Dental 510(k) Summary

May 06, 2011

Sirona Dental CAD/CAM System

APPENDIX H • Page 1 of 32

#### **4 Intended Use**

The Sirona Dental CAD/CAM System is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations. For the titanium bases SSO 3.5 L and SBL 3.3 L, the indication is restricted for replacement of single lateral incisors in the maxilla and lateral and central incisors in the mandible. The system consists of three major parts: TiBase, inCoris mesostructure, and CAD/CAM software. Specifically, the inCoris mesostructure and TiBase components make up a two-piece abutment which is used in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity. The inCoris mesostructure may also be used in conjunction with the Camlog Titanium base CAD/CAM (types K2244.xxxx) (K083496) in the Camlog Implant System. The CAD/CAM software is intended to design and fabricate the inCoris mesostructure. The inCoris mesostructure and TiBase two-piece abutment is compatible with the following implants systems:

- Nobel Biocare Replace (K020646)
- Nobel Biocare Branemark (K022562)
- Friadent Xive (K013867)
- Biomet 3i Osseotite (K980549)
- Astra Tech Osseospeed (K091239)
- Zimmer Tapered Screw-Vent (K061410)
- Straumann SynOcta (K061176)
- Straumann Bone Level (K053088, K062129, K060958)
- Biomet 3i Certain (K014235, K061629)
- Nobel Biocare Active (K071370)

#### **5 Device Description**

The Sirona Dental CAD/CAM-System takes optical impressions and records the topographical characteristics of teeth, dental impressions, or stone models. Dental restorative prosthetic devices are manufactured using computer aided design and fabrication. The system also features the processing of mesostructures, a dental restorative prosthetic device used in conjunction with endosseous dental implant abutments.

The system that features the processing of mesostructures comprises

- Titanium bases TiBase and Camlog

- inCoris ZI meso blocks
- Sirona Dental CAD/CAM Design and fabricating devices

Titanium bases are used as an implant prosthetic titanium base for adhesion to mesostructures to restore function and aesthetics in the oral cavity.

inCoris ZI meso blocks are used in manufacturing individually designed mesostructures, which are glued to a fitting titanium base after milling and sintering.

Sirona Dental CAD/CAM design and fabricating devices feature the processing of mesostructures, a dental restorative prosthetic device used in conjunction with endosseous dental implant abutments, i.e. it is an accessory to it. This component consists of the devices CEREC3, CEREC AC, inEos, inEos Blue, CEREC MCXL and inLab MCXL.

## **5.1 TiBase**

### **5.1.1 Device Function**

The Sirona TiBase is a premanufactured prosthetic component directly connected to dedicated endosseous dental implants and is intended for use as an aid in prosthetic rehabilitation.

The Sirona offering consists of the titanium base TiBase, the abutment screw and the scanbody. The parts are marketed non-sterile and for single use only.

The Sirona TiBase is bonded to an individually designed mesostructure, a ceramic prosthetic/restoration, that supports the final restoration. The mesostructure is milled from an inCoris ZI meso block with Sirona CAD/CAM milling machines CEREC or inLab, and sintered afterwards.

The two piece abutment is mounted onto the implant and fixed with a screw.

The scope of delivery contains a scanbody (ABS plastic) which is mounted on a TiBase in order to acquire the topographical surface of the area where the endosseous dental implant abutment is located with Sirona Dental CAD/CAM fabricating devices. From the acquired data the position of the implant can be calculated. After an optical impression has been taken the scanbody is removed.

Sirona TiBase devices are compatible with following systems (Table 1):

**Table 1: Sirona TiBase Devices Compatibility**

Sirona TiBase	Compatible System		
	Manufacturer	System	Diameter
NBRS 3.5	Nobel Biocare	Replace® NP	3,5 mm
NBRS 4.3		Replace® RP	4.3 mm
NBRS 5.0		Replace® WP	5.0 mm
NBRS 6.0		Replace® 6.0	6.0 mm
NBB 3.4	Nobel Biocare	Brånemark®	3.4 mm
NBB 4.1		Brånemark®	4.1 mm
SSO 3.5	Straumann	Tissue level NN	3.5 mm
SSO 4.8		Tissue level RN	4.8 mm
SSO 6.5		Tissue level WN	6.5 mm
ATOS 3.5/4.0	Astra Tech	OsseoSpeed™	3.5 S / 4.0 S mm
ATOS 4.5/5.0		OsseoSpeed™	4.5 / 5.0 mm
FX 3.4	Friadent	Frialit® / Xive®	3.4 mm
FX 3.8		Frialit® / Xive®	3.8 mm
FX 4.5		Frialit® / Xive®	4.5 mm
FX 5.5		Frialit® / Xive®	5.5 mm
BO 3.4	Biomet 3i	Osseotite (Connec-tion type: Ex. Hex)	3.4 mm
BO 4.1		Osseotite (Connec-tion type: Ex. Hex)	4.1 mm

Sirona TiBase	Compatible System		
	Manufacturer	System	Diameter
BO 5.0		Osseotite (Connection type: Ex. Hex)	5.0 mm
ZTSV 3.5	Zimmer	Tapered Screw-Vent®	3.5 mm
ZTSV 4.5		Tapered Screw-Vent®	4.5 mm
ZTSV 5.7		Tapered Screw-Vent®	5.7 mm
NB A 4.5	Nobel Biocare	Nobel Active NP	3.5mm
NB A 5.0		Nobel Active NP	4.3 / 5.0mm
S BL 3.3	Straumann®	Bone Level NC	3.3mm
S BL 4.1		Bone Level NC	4.1 / 4.8mm
B C 3.4	Biomet 3i	Certain®	3.4mm
B C 4.1		Certain®	4.1mm
B C 5.0		Certain®	5.0mm

### 5.1.2 Scientific Concept

The underlying scientific concept is the use of an already introduced technology of a titanium base abutment combined with individually CAD/CAM fabricated ceramic prosthetics.

### 5.1.3 Physical and Performance Characteristics

#### 5.1.3.1 Design

The TiBase devices have various diameters, are compatible with dedicated implant systems, and fit to compatible implants as provided in Table 2.



**Table 2: Implant Compatibility**

<b>TiBase</b>	<b>Implant- Manufacturer</b>	<b>Implant-System</b>	<b>510(k) Implant</b>
NBRS	Nobel Biocare	Replace	K020646
NBB	Nobel Biocare	Branemark	K022562
FX	Friadent	Xive	K013867
BO	Biomet 3i	Osseotite	K980549
ATOS	Astra Tech	OsseoSpeed	K091239
ZTSV	Zimmer	Tapered Screw-Vent	K061410
SSO	Straumann	SynOcta	K061176
NB A	Nobel Biocare	Nobel Active	K071370
S BL	Straumann®	Bone Level	K053088
B C	Biomet 3i	Certain®	K014235

**5.1.3.2 Material Used**

TiBase and abutment screw are made of Ti6Al4V.

**5.1.3.3 Physical Properties**

TiBase material composition and mechanical properties comply with ISO 5832-3:1996, Implants for surgery -- Metallic materials -- Part 3: Wrought titanium 6-aluminium 4-vanadium alloy

Sirona TiBase devices are compatible with systems listed in Table 1.

**5.2 inCoris ZI meso**

**5.2.1 Device Description**

The inCoris ZI meso offerings are blocks of various sizes from which individual dental mesostructures are grinded by milling machines (inLab MCXL, CEREC MCXL). The mesostructure is a part of a 2 part endosseous dental implant

abutment which comprises a titanium base and a zirconium oxide mesostructure. The connection geometries are prefabricated.

### 5.2.2 Scientific Concept

The underlying scientific concept is the use of an already introduced technology of a titanium base abutment combined with individually CAD/CAM fabricated ceramic prosthetics.

### 5.2.3 Physical and Performance Characteristics

#### 5.2.3.1 Design

The inCoris ZI meso are blocks of various sizes. The marketed ceramic is pre-sintered. One side of a block is mounted to a mandrel that will be inserted in the spindle's clamping chuck of the grinding machine. The connection geometry to titanium bases is prefabricated, i.e. already included in the shipped block. Connection geometries fit on Camlog (type K2244.xxxx) and Sirona (Tibase) titanium bases (Table 3 and Table 4).

**Table 3: Sirona inCoris ZI meso - TiBase Devices Compatibility**

Titanium Base			Ceramic Block		
TiBase (Sirona)	REF	Diameter	inCoris ZI meso (Sirona)	REF	Color inCoris ZI meso (Sirona)
NBRS 3.5	6282474	3,5 mm	inCoris ZI meso L	62 31 810	F0.5
				62 31 836	F2
NBRS 4.3	6282482	4.3 mm	inCoris ZI meso L	62 31 810	F0.5
				62 31 836	F2
NBRS 5.0	6282490	5.0 mm	inCoris ZI meso L	62 31 810	F0.5
				62 31 836	F2
NBRS 6.0	6282508	6.0 mm	inCoris ZI meso L	62 31 810	F0.5
				62 31 836	F2
NBB 3.4	6282516	3.4 mm	inCoris ZI	62 31 810	F0.5

Titanium Base			Ceramic Block		
TiBase (Sirona)	REF	Diameter	inCoris ZI meso (Sirona)	REF	Color inCoris ZI meso (Sirona)
			meso L	62 31 836	F2
NBB 4.1	6282524	4.1 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
SSO 3.5	6284231	3.5 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
SSO 4.8	6284249	4.8 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
SSO 6.5	6284256	6.5 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
ATOS 3.5/4.0	6282532	3.5 S / 4.0 S mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
ATOS 4.5/5.0	6282540	4.5 / 5.0 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
FX 3.4	6282433	3.4 mm	inCoris ZI meso S	62 31 802 62 31 828	F0.5 F2
FX 3.8	6282441	3.8 mm	inCoris ZI meso S	62 31 802 62 31 828	F0.5 F2
FX 4.5	6282458	4.5 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
FX 5.5	6282466	5.5 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
BO 3.4	6282557	3.4 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2

Titanium Base			Ceramic Block		
TiBase (Sirona)	REF	Diameter	inCoris ZI meso (Sirona)	REF	Color inCoris ZI meso (Sirona)
BO 4.1	6282565	4.1 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
BO 5.0	6282573	5.0 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
ZTSV 3.5	6282581	3.5 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
ZTSV 4.5	6282599	4.5 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
ZTSV 5.7	6282607	5.7 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
NB A 4.5	6308188	3.5mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
NB A 5.0	6308253	4.3 / 5.0mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
S BL 3.3	6308154	3.3mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
S BL 4.1	6308337	4.1 / 4.8mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
B C 3.4	6308048	3.4mm	inCoris ZI meso S	62 31 802 62 31 828	F0.5 F2
B C 4.1	6308097	4.1mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2

Titanium Base			Ceramic Block		
TiBase (Sirona)	REF	Diameter	inCoris ZI meso (Sirona)	REF	Color inCoris ZI meso (Sirona)
B C 5.0	6308121	5.0mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2

**Table 4: Sirona inCoris ZI meso - Camlog Devices Compatibility**

Titanium Base			Ceramic Block		
Camlog	REF	Dia- meter	inCoris ZI meso (Sirona)	REF	Color inCoris ZI meso (Sirona)
K2244.3348	K2244.3348	3.3	inCoris ZI meso S	62 31 802 62 31 828	F0.5 F2
K2244.3848	K2244.3848	3.8	inCoris ZI meso S	62 31 802 62 31 828	F0.5 F2
K2244.4348	K2244.4348	4.3	inCoris ZI meso S	62 31 802 62 31 828	F0.5 F2
K2244.5048	K2244.5048	5.0	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
K2244.6048	K2244.6048	6.0	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2

### 5.2.3.2 Material Used

The inCoris ZI meso are pre-sintered zirconium oxide ceramic blocks. The metal block holder is made of aluminum. The material is composed of:

ZrO <sub>2</sub> +HfO <sub>2</sub> +Y <sub>2</sub> O <sub>3</sub>	> 99.0%
Al <sub>2</sub> O <sub>3</sub>	< 0.5%
Other oxides	< 0.5%

### 5.2.3.3 Physical Properties

The final technical data of inCoris ZI meso are (after final sintering):

Density:	6.06 g/cm <sup>3</sup>
Coefficient of thermal expansion (CTE):	11.0*10 <sup>-6</sup> K <sup>-1</sup>
Flexural strength:	> 900 MPa
Fracture toughness (KIC):	5.9 MPa·m <sup>1/2</sup>

## 5.3 Sirona Dental CAD/CAM Design and fabrication Devices

### 5.3.1 Device Description

The Sirona Dental CAD/CAM Design and fabricating devices for processing mesostructures includes

- Optical acquisition or recording of the topographical characteristics of dental impressions, or stone models using the devices Acquisition unit CEREC 3, CEREC AC, and stationary scanning system inEos Blue
- Design of mesostructures and processing the acquired or recorded data for these purposes using Sirona Dental CAD/CAM Software which runs on a CEREC 3, CEREC AC or PC. Design is performed by a dentist or dental technician
- Milling of the mesostructure using CEREC MCXL or inLab MCXL milling machines from ceramic blocks intended for dental restorations and mesostructures

The Sirona Dental CAD/CAM Design and fabricating devices also processes other dental restorations like crowns, bridge-frameworks, inlays, onlays all regulated under 21 CFR 872.3661, Optical Impression Systems for CAD/CAM, for such intended use.

### **5.3.2 Scientific Concept**

The underlying scientific concept is the use of CAD/CAM technology for the optical acquisition of the topographical characteristics of dental impressions, and models, the design of individual mesostructures using recorded data (CAD), and eventually fabricating (milling) these designed mesostructures (CAM).

### **5.3.3 Physical and Performance Characteristics**

#### **5.3.3.1 Design**

Acquisition unit: the device consists of a camera for acquiring optical topographical characteristics of dental impressions. The recorded data are used for the design of individual mesostructures using CAD techniques specific to the dental field.

Fabricating devices: the devices mill the individual designed mesostructures from incoris ZI meso blocks. For this purpose the chucked block and the milling tools move according to prescribed trajectories to generate the shape which is intended to be milled.

#### **5.3.3.2 Materials Used**

Not applicable.

#### **5.3.3.3 Physical Properties**

Not applicable.

## **6 Summary of the technological characteristics**

### **6.1 TiBase**

All proposed and predicate titanium bases and screws are made of Ti6Al4V, medical grade 5. Connection interfaces to the implants are identical for each defined diameter and connection type. Connection interfaces to dental restorations differ in that proposed devices have an additional notch.

An extensive list is provided in Table 5.

### **6.2 inCoris ZI meso**

Proposed and predicate (K100152) device are inCoris ZI meso.

There has been no modification to inCoris ZI meso from Premarket Notification K100152. Specifically, the following aspects remain identical:

- Composition
- Material properties
- Thickness / design restrictions
- Shape and bonding surface of connection interfaces to Camlog and Sirona Tibase for Sirona inCoris ZI meso blocks
- Bonding material

InCoris ZI meso is bonded to titanium bases for supporting further dental restorations.

InCoris ZI meso material is made of zirconium oxide. The composition of inCoris complies with ISO standard 13356:1997, "Implants for surgery, Ceramic materials based on yttria-stabilized tetragonal zirconia (Y-TZP)". Such compositions are classified as biocompatible by this standard.

An extensive list is provided in Table 6.

### **6.3 Sirona Dental CAD/CAM Design and fabrication Devices**

There has been no modification to the Sirona Dental CAD/CAM Design and fabricating devices from Premarket Notification K100152. Specifically, the following aspects remain identical:

- optical impressions record topographical characteristics of teeth, dental impressions, or stone models for use in the computer aided design and fabrication of dental restorative prosthetic devices in conjunction with endosseous dental implant abutments, i.e. it is an accessory
- features the transfer of data of the optical impression to a remote milling machine via internet or exportation/importation of milling data

The software database of titanium bases has been extended to cover new additional titanium bases.

An extensive list is provided in Table 7.



**Table 5: Comparison of Sirona TiBase to Predicate Devices**

Proposed Device	Predicate Device						Screw geometry			
	Manufacturer	System	Diameter	Titanium Base	Predicate Devices K-Number	Abutment and Screw made of Ti6Al4V		Identical connection geometry to abutments	Connection geometry	Anti-rotational features
TiBase										
NBRS 3.5	Nobel Biocare	Replace® NP	3.5 mm	Nobel Biocare product catalog page 14. Product-No. 32376	K091756	yes	yes	Internal 3 te- nons	yes	same
NBRS 4.3	Nobel Biocare	Replace® RP	4.3 mm	Product-No. 32377	K091756	yes	yes	Internal 3 te- nons	yes	same

May 06, 2011

Sirona Dental 510(k) Summary

Sirona Dental CAD/CAM System

APPENDIX H • Page 14 of 32

Proposed Device	Predicate Device						Identical connection geometry to abutments	Connection geometry Type	Connection geometry	Screw geometry
	Manufacturer	System	Diameter	Titanium Base	Predicate Devices K-Number	Abutment and Screw made of Ti6Al4V				
TiBase										
NBRS 5.0	Nobel Biocare	Replace® WP	5.0 mm	Product-No. 32378	K091756	yes	Internal 3 te-nons	yes	same	
NBRS 6.0	Nobel Biocare	Replace® 6.0	6.0 mm	Product-No. 32375	K091756	yes	Internal 3 te-nons	yes	same	

May 06, 2011

Sirona Dental 510(k) Summary

Sirona Dental CAD/CAM System

Proposed Device TIBase	Predicate Device						Abutment and Screw made of Ti6Al4V	Identical connection geometry to abutments	Connection geometry Type	Connection geometry	Screw geometry
	Manufacturer	System	Diameter	Titanium Base	Predicate Devices K-Number						
NBB 3.4	Nobel Biocare	Brånemark®	3.4 mm	Nobel Biocare product catalog page 14. Product-No. 32396	K091756		yes	External Hexagonal	yes	same	
NBB 4.1	Nobel Biocare	Brånemark®	4.1 mm	Product-No. 32397	K091756		yes	External Hexagonal	yes	same	

May 06, 2011

Sirona Dental 510(k) Summary

Sirona Dental CAD/CAM System

APPENDIX H • Page 16 of 32

Proposed Device	Predicate Device					Identical connection geometry to abutments	Connection geometry Type	Connection geometry	Screw geometry
	Manufacturer	System	Diameter	Titanium Base	Predicate Devices K-Number				
TIBase						Abutment and connection geometry of Ti6Al4 V			
SSO 3.5	Straumann	Tissue level NN	3.5 mm	Straumann product catalog page 51. Product-No. 048.505	K081005	yes	external octagonal	yes	same
SSO 4.8	Straumann	Tissue level RN	4.8 mm	Product-No. 048.600 (p 55)	K081005	yes	Internal Octagonal	yes	same
SSO 6.5	Straumann	Tissue level WN	6.5 mm	Product-No. 048.606	K081005	yes	Inter-	yes	same

May 06, 2011

Sirona Dental 510(k) Summary

Sirona Dental CAD/CAM System

Proposed Device TIBase	Predicate Device					Abutment and Screw made of Ti6Al4V	Identical connection geometry to abutments	Connection geometry Type	Connection geometry Anti-rotational features	Screw geometry
	Manufacturer	System	Diameter	Titanium Base	Predicate Devices K-Number					
					(p 62)					
ATOS 3.5/4.0	Astra Tech	OsseoSpeed™	3.5 S / 4.0 S mm	Product-No. 24285	K081666	yes	Internal Hexagonal	yes	same	
ATOS 4.5/5.0	Astra Tech	OsseoSpeed™	4.5 / 5.0 mm	Product-No. 24235	K081666	yes	Internal Hexagonal	yes	same	

May 06, 2011

Sirona Dental 510(k) Summary  
Sirona Dental CAD/CAM System

Proposed Device	Predicate Device						Abutment and connection geometry of Ti6Al4V	Identical connection geometry to abutments	Connection geometry Type	Connection geometry	Screw geometry
	Manufacturer	System	Diameter	Titanium Base	Predicate Devices K-Number	Abutment and connection geometry of Ti6Al4V					
TI Base											
FX 3.4	Friadent	Frialit® / Xive®	3.4 mm	Friadent product catalog page 32. Product-No. 46-2132	K032158	yes	yes	Internal Hexagonal	yes	same	
FX 3.8	Friadent	Frialit® / Xive®	3.8 mm	Product-No. 46-2142	K032158	yes	yes	Internal Hexagonal	yes	same	

May 06, 2011

Sirona Dental 510(k) Summary

Sirona Dental CAD/CAM System

Proposed Device	Predicate Device					Identical connection geometry to abutments	Connection geometry Type	Connection geometry	Screw geometry
	TI Base	Manufacturer	System	Diameter	Titanium Base				
FX 4.5	Friadent	Frialit® / Xive®	4.5 mm	Product-No. 46-2152	K032158	yes	Internal Hexagonal	yes	same
FX 5.5	Friadent	Frialit® / Xive®	5.5 mm	Product-No. 46-2162	K032158	yes	Internal Hexagonal	yes	same
BO 3.4	Biomet 3i	Osseotite (Connection)	3.4 mm	3i Biomet product catalog page 25.	K072642	yes	External	yes	same

Sirona Dental 510(k) Summary

May 06, 2011

Sirona Dental CAD/CAM System

APPENDIX H • Page 20 of 32

Proposed Device	Predicate Device						Abutment and Screw made of Ti6Al4V	Identical connection geometry to abutments	Connection geometry Type	Connection geometry	Screw geometry
	Manufacturer	System	Diameter	Titanium Base	Predicate Devices K-Number						
TiBase											
BO 4.1	Biomet 3i	Osseointe (Connec-tion type: Ex. Hex)	4.1 mm	Product-No. APP452G	K072642	Product-No. MAP32G	yes	yes	Hexagonal External Hexagonal	yes	same
BO 5.0	Biomet 3i	Osseointe (Connec-tion type: Ex. Hex)	5.0 mm	Product-No. WPP552G	K072642		yes	yes	External Hexagonal	yes	same

Sirona Dental 510(k) Summary

Sirona Dental CAD/CAM System

May 06, 2011

APPENDIX H • Page 21 of 32



Proposed Device TIBase	Predicate Device				Abutment and Screw made of Ti6Al4V	Identical connection geometry to abutments	Connection geometry Type	Connection geometry Anti-rotational features	Screw geometry
	Manufacturer	System	Diameter	Titanium Base Titanium Base Zimmer product catalog page 9. Product-No. ZOA342S Product-No. ZOA442S					
ZTSV 3.5	Zimmer	Tapered Screw-Vent®	3.5 mm	Zimmer product catalog page 9. Product-No. ZOA342S	K060880	yes	Internal Hexagonal	yes	same
ZTSV 4.5	Zimmer	Tapered Screw-Vent®	4.5 mm	Product-No. ZOA442S	K060880	yes	Internal Hexagonal	yes	same

May 06, 2011

Sirona Dental 510(k) Summary  
Sirona Dental CAD/CAM System

Proposed Device TiBase	Predicate Device						Abutment and Screw made of Ti6Al4V	Identical connection geometry to abutments	Connection geometry Type	Connection geometry Anti-rotational features	Screw geometry
	Manufacturer	System	Diameter	Titanium Base	Predicate Devices K-Number						
ZTSV 5.7	Zimmer	Tapered Screw-Vent®	5.7 mm	Product-No. ZOAS62S	K060880	yes	yes	Internal Hexagonal	yes	same	
NB A 4.5	Nobel Biocare	Nobel Active NP	3.5mm	Nobel Biocare product Catalog page 71. Product-No. 34194	K102436	yes	yes	Internal Hexagonal	yes	same	

May 06, 2011

Sirona Dental 510(k) Summary

Sirona Dental CAD/CAM System

Proposed Device TiBase	Predicate Device					Abutment and Screw made of Ti6Al4V	Identical connection geometry to abutments	Connection geometry Type	Connection geometry	Screw geometry
	Manufacturer	System	Diameter	Titanium Base	Predicate Devices K-Number					
NB A 5.0	Nobel Biocare	Nobel Active RP	4.3 / 5.0mm	Product-No. 34198	K102436	yes	yes	Internal Hexagonal	yes	same
S BL 3.3 ®	Straumann	Bone Level NC	3.3mm	Straumann product Catalog Product-No. 022.2102	K062129	yes	yes	Internal 4 slots	yes	same

May 06, 2011

Sirona Dental 510(k) Summary

Sirona Dental CAD/CAM System

Proposed Device TIBase	Predicate Device					Abutment and Screw made of Ti6Al4V	Identical connection geometry to abutments	Connection geometry Type	Connection geometry	Screw geometry
	Manufacturer	System	Diameter	Titanium Base	Predicate Devices K-Number					
S BL 4.1	Straumann	Bone Level RC	4,1 / 4,8mm	Product-No. 022.4102	K062129	yes	yes	Internal 4 slots	yes	same
B C 3.4	Biomet 3i	Certain®	3,4mm	Biomet product Catalog page 5 Product-No. IMAP32G	K073345	yes	yes	Internal Hexagonal	yes	same
B C 4.1	Biomet 3i	Certain®	4,1 mm	Biomet product	K073345	yes	yes	Internal	yes	same

May 06, 2011

Sirona Dental 510(k) Summary

Sirona Dental CAD/CAM System

Proposed Device TIBase	Predicate Device					Abutment and Screw made of Ti6Al4V	Identical connection geometry to abutments	Connection geometry Type	Connection geometry	Screw geometry
	Manufacturer	System	Diameter	Titanium Base	Predicate Devices K-Number					
				Catalog page 5 Product-No. IAPP452G				Hexagonal		
B C 5.0	Biomet 3i	Certain®	5.0mm	Biomet product Catalog page 5 Product-No. IWPP562G	K073345	yes	Internal Hexagonal	yes	same	

May 06, 2011

Sirona Dental 510(k) Summary

Sirona Dental CAD/CAM System

APPENDIX H • Page 26 of 32

**Table 6: Comparison of InCoris ZI meso to Predicate Device**

	InCoris ZI meso	InCoris ZI meso (K100152)
Intended use	inCoris ZI meso blocks are used in manufacturing individually designed mesostructures, which are glued to a fitting titanium base after milling and sintering.	inCoris ZI meso blocks are used in manufacturing individually designed mesostructures, which are glued to a fitting titanium base after milling and sintering.
Application	inCoris ZI mesostructures can only be used for the intended titanium bases or implants. Allocation of the connection size to the respective titanium base can be determined by the scanbody set of the respective implant system.  Please observe the indications and contraindications of the implant.	inCoris ZI mesostructures can only be used for the intended titanium bases or implants. Allocation of the connection size to the respective titanium base can be determined by the scanbody set of the respective implant system.  Please observe the indications and contraindications of the implant.
Contra-Indications	<ul style="list-style-type: none"> <li>• Insufficient oral hygiene</li> <li>• Insufficient space available</li> <li>• Bruxism</li> <li>• For mesostructure-geometry with angulation correction greater than 20° to the implant axis</li> <li>• For mesostructure-geometry with angulation correction to the implant axis for Camlog only</li> <li>• For individual tooth restorations with free end saddle</li> <li>• For restorations with a length to implant length ratio of more than 1:1.25</li> </ul>	<ul style="list-style-type: none"> <li>• Insufficient oral hygiene</li> <li>• Insufficient space available</li> <li>• Bruxism</li> <li>• For restorations with angulation correction to the implant axis</li> <li>• For individual tooth restorations with free end saddle</li> <li>• For restorations with a length to implant length ratio of more than 1:1.25</li> </ul>
Technical Data		
Block-Material Composition	<ul style="list-style-type: none"> <li>• ZrO<sub>2</sub>+HfO<sub>2</sub>+Y<sub>2</sub>O<sub>3</sub>: &gt; 99.0%</li> <li>• Y<sub>2</sub>O<sub>3</sub>: 5.2%</li> <li>• HfO<sub>2</sub>: 2%</li> </ul>	<ul style="list-style-type: none"> <li>• ZrO<sub>2</sub>+HfO<sub>2</sub>+Y<sub>2</sub>O<sub>3</sub>: &gt; 99.0%</li> <li>• Y<sub>2</sub>O<sub>3</sub>: 5.2%</li> <li>• HfO<sub>2</sub>: 2%</li> </ul>

	InCoris ZI meso	InCoris ZI meso (K100152)
	<ul style="list-style-type: none"> <li>Al<sub>2</sub>O<sub>3</sub>: ≤ 0.05%</li> <li>Fe<sub>2</sub>O<sub>3</sub>: 0.3%</li> </ul>	<ul style="list-style-type: none"> <li>Al<sub>2</sub>O<sub>3</sub>: ≤ 0.05%</li> <li>Fe<sub>2</sub>O<sub>3</sub>: 0.3%</li> </ul>
Density (sintered)	6.06 g/cm <sup>3</sup>	6.06 g/cm <sup>3</sup>
Coefficient of thermal expansion (CTE)	11.0*10 <sup>-6</sup> K <sup>-1</sup> (20 °C - 500 °C)	11.0*10 <sup>-6</sup> K <sup>-1</sup> (20 °C - 500 °C)
Flexural strength	> 900MPa	> 900MPa
Fracture toughness (K <sub>IC</sub> )	5.8 MPa·m <sup>1/2</sup>	5.8 MPa·m <sup>1/2</sup>
Grain Size	about 0.5 μm	about 0.5 μm
Bonding Material	Panavia F 2.0 (www.kuraray-dental.de)	Panavia F 2.0 (www.kuraray-dental.de)

**Table 7: Comparison of Sirona Dental CAD/CAM Design and fabricating Devices to Predicate Devices**

	<b>Sirona Dental CAD/CAM Design and fabricating Devices</b>	<b>Sirona Dental CAD/CAM Design and fabricating Devices (K100152)</b>
Used for	<p>The Sirona Dental CAD/CAM-System is indicated for taking optical impressions to record the topographical characteristics of teeth, dental impressions, or stone models by computer aided design and fabrication in patients that require dental restorative prosthetic devices. The system also features the processing of mesostructures, a dental restorative prosthetic device used in conjunction with endosseous dental implant abutments.</p> <p>Sirona Dental CAD/CAM Design and fabricating devices feature the processing of mesostructures, a dental restorative prosthetic device used in conjunction with endosseous dental implant abutments, i.e. it is an accessory to it. Devices which feature the processing of mesostructures comprises CEREC3, CEREC AC, inEos, inEos Blue, CEREC MCXL and inLab MCXL</p>	<p>The Sirona Dental CAD/CAM-System is indicated for taking optical impressions to record the topographical characteristics of teeth, dental impressions, or stone models by computer aided design and fabrication in patients that require dental restorative prosthetic devices. The system also features the processing of mesostructures, a dental restorative prosthetic device used in conjunction with endosseous dental implant abutments.</p> <p>Sirona Dental CAD/CAM Design and fabricating devices feature the processing of mesostructures, a dental restorative prosthetic device used in conjunction with endosseous dental implant abutments, i.e. it is an accessory to it. Devices which feature the processing of mesostructures comprises CEREC3, CEREC AC, inEos, inEos Blue, CEREC MCXL and inLab MCXL</p>
Used with	Sirona Dental CAI/CAM Hardware	Sirona Dental CAI/CAM Hardware
Controlling of recording process (CAI)  (optical impression)	Yes	Yes
Processing the recorded	Yes	Yes



	Sirona Dental CAD/CAM Design and Fabricating Devices	Sirona Dental CAD/CAM Design and Fabricating Devices (1300152)
data (data of optical impression) (CAD)		
Export of milling data to milling machine	Yes	Yes
Administration of patient data	Yes	Yes
Further functions	Calibration of CAI/CAM hardware	Calibration of CAI/CAM hardware
Online capability	Option to upload/download the data from a web portal (Cerec Connect), to have CAI and CAM operating on two different locations connected via Internet	Option to upload/download the data from a web portal (Cerec Connect), to have CAI and CAM operating on two different locations connected via Internet
Scan Implant Interface/surface	Yes  (or with mounted scanbody)	Yes  (or with mounted scanbody)
Scan custom wax-up	Yes	Yes
Preparation of individual restoration ("meso-structure") to be mounted on the abutment	Yes	Yes
Bond of milled zirconia/ceramic individual meso-structure to metal abutment	Yes	Yes
Create of fitting crown to be mounted on top of meso-structure	Yes	Yes
Used with	Sirona Dental CAI/CAM Hardware	Sirona Dental CAI/CAM Hardware

	Sirona Dental CAD/CAM Design and Fabricating Devices	Sirona Dental CAD/CAM Design and Fabricating Devices (13100152)
Used for	CAD creation of dental restorations including inlays, onlay, veneers, crowns, bridges and meso-structure to be mounted on top of abutments	CAD creation of dental restorations including inlays, onlay, veneers, crowns, bridges and meso-structure to be mounted on top of abutments
Controlling of measurement process (CAI)	Yes	Yes
Processing the measurement data (CAD)	Yes	Yes
Export to milling machine	Yes	Yes
Administration of patient data	Yes	Yes
Further functions	Calibration of CAI/CAM hardware	Calibration of CAI/CAM hardware

## 7 Nonclinical Testing

According to FDA Guidance "Class II Special Controls Guidance Document: Root-form Endosseous Dental Implants and Endosseous Dental Implant Abutments", May 12, 2004, fatigue testing has been performed for angled abutments. A reverse-engineering analysis has demonstrated that the devices are identical with their predicates and compatible with their mating implants.

Software validation testing has been performed according to IEC 62304:2006. A warning has been added warning the user that abutments with an angle of greater than 20° are out of specification.

## 8 Clinical Testing

Clinical testing is not required and has not been performed.

## **9 Conclusion**

Based on a comparison of intended use, indications, construction materials, principal of operations, features and technical data, the Sirona Dental CAD/CAM System which comprises of titanium bases TiBase, inCoris ZI meso blocks and Sirona Dental CAD/CAM Design and fabricating devices are safe and effective their intended use and perform as well as and are substantially equivalent to their Predicate Devices.

---

Sirona Dental 510(k) Summary      May 06, 2011

Sirona Dental CAD/CAM System

APPENDIX H • Page 32 of 32



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - W066-G609  
Silver Spring, MD 20993-0002

Mr. Fritz Kolle  
Sirona Dental Systems GmbH  
Fabrikstrasse 31  
Bensheim  
Germany D-64625

FEB 17 2012

Re: K111421  
Trade/Device Name: Sirona Dental CAD/CAM System  
Regulation Number: 21 CFR 872.3630  
Regulation Name: Endosseous Dental Implant Abutment  
Regulatory Class: II  
Product Code: NHA  
Dated: February 14, 2012  
Received: February 16, 2012

Dear Mr. Kolle:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Page 2- Mr. Kollé

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

  
Anthony D. Watson, BS, MS, MBA  
Director  
Division of Anesthesiology, General Hospital,  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number (if known):

Device Name: Sirona Dental CAD/CAM System

Indications for Use:

The Sirona Dental CAD/CAM System is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations. The system consists of three major parts: TiBase, InCoris mesostructure, and CAD/CAM software. Specifically, the InCoris mesostructure and TiBase components make up a two-piece abutment which is used in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity. The InCoris mesostructure may also be used in conjunction with the Camlog Titanium base CAD/CAM (types K2244.xxxx) (K083496) in the Camlog Implant System. The CAD/CAM software is intended to design and fabricate the InCoris mesostructure. The InCoris mesostructure and TiBase two-piece abutment is compatible with the following implant systems:

- Nobel Biocare Replace (K020646)
- Nobel Biocare Branemark (K022562)
- Friadent Xive (K013867)
- Biomet 3i Osseotite (K980549)
- Astra Tech Osseospeed (K091239)
- Zimmer Tapered Screw-Vent (K061410)
- Straumann SynOcta (K061176)
- Straumann Bone Level (K053088)
- Biomet 3i Certain (K014235)
- Nobel Biocare Active (K071370)

*Y. S. Batzmas for Dr. S. Runner*  
(Division Sign-Off)  
Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K111421

Prescription Use X  
(Part 21 CFR 801 Subpart D)

OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NECESSARY)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Sirona Dental Systems 510(k)  
Sirona Dental CAD/CAM System

May 06, 2011

Page x

**Traditional 510(k) Submission**  
**Straumann® Variobase™ Abutments**

Appendix 23

---

**Appendix 23 – Lava™ Frame, Lava™ Frame Shade 510(k)  
Summary**

**510(k) SUMMARY OF SAFETY AND EFFECTIVENESS**

**Submitter**

Company: .....3M ESPE AG  
Street: .....ESPE Platz  
ZIP-Code, City:.....D-82229 Seefeld  
Federal State: .....Bavaria  
Country: .....Germany  
Establishment Registration Number .....9611385  
Official Correspondent: .....Dr. Andreas Petermann,  
.....Manager Regulatory Affairs  
Phone: .....011-49-8152-700 1395  
Fax: .....011-49-8152-700 1869  
E-mail:.....Andreas.Petermann@mmm.com  
Date:.....July 23, 2007

**Name of Device**

Proprietary Name:.....Lava™ Frame, Lava™ Frame Shade  
Classification Name:.....Porcelain powder for clinical use  
.....Endosseous dental implant abutment  
Common Name: .....All-ceramic core material  
.....All-ceramic stain solution  
.....Abutment

**Predicate Device**

ALTATEC Camlog Implant System and  
Abutments by Altatec Biotechnologies .....K032448



K072055

11562

Description for the Premarket Notification

Lava™ abutment made from Lava™ Frame zirconia mill blanks and dyed with Lava™ Frame Shade is classified as endosseous dental implant abutment (21 C.F.R. § 872.3630) because it is a prosthetic component directly connected to the endosseous dental implant and is intended for use as an aid in prosthetic rehabilitation.

Lava™ Frame and Lava™ Frame Shade are parts of the Lava™ system (K011394). Lava™ Frame Zirconia mill blanks are used for the fabrication of frameworks for all-ceramic restorations. The frameworks for onlays, inlays, veneers, crowns and bridges are designed and manufactured by CAD/CAM technology, whereas the CAM fabricated Lava™ Abutments made from Lava™ Frame Zirconia mill blanks will be designed by means of a traditional wax up technique. The wax up will be scanned (Lava™ Scan, K062493) and milled without any further design step in the CNC milling unit Lava™ Form. After milling, the abutments are dyed with one of the 7 Lava™ Frame Shade dyeing liquids as required to achieve the desired tooth color, then sintered. The dyed abutments are sintered using the specialized program of the Lava™ Therm sintering furnace.

The wax up designed abutment will be cemented to a titanium interface which will be subsequently screwed into the respective implant (e.g. Camlog, Altatec Biotechnologies).

The comparison for composition, performance data and indications for use shows that Lava™ abutment made from Lava™ Frame and dyed with Lava™ Frame Shade is substantially equivalent to the predicate device.

In summary, it can be concluded that safety and effectiveness requirements for Lava™ Frame and Lava™ Frame Shade for the fabrication of abutments are completely met.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room –WO66-G609  
Silver Spring, MD 20993-0002

Dr. Desi W. Soegiarto  
Regulatory Affairs Specialist  
3M ESPE AG Dental Products  
ESPE Platz  
Seefeld, Bavaria  
GERMANY D-82229

MAR 29 2011

Re: K072055  
Trade/Device Name: Lava™ Frame, Lava™ Frame Shade  
Regulation Number: 21 CFR 872.3630  
Regulation Name: Endosseous Dental Implant Abutment  
Regulatory Class: II  
Product Code: NHA  
Dated: February 18, 2008  
Received: February 21, 2008

Dear Dr. Soegiarto:

This letter corrects our substantially equivalent letter of February 26, 2008.

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Page 2– Dr. Soegiarto

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Anthony D. Watson". To the left of the signature is a small, stylized handwritten mark that looks like "fa".

Anthony D. Watson, B.S., M.S., M.B.A.  
Director

Division of Anesthesiology, General Hospital,  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K072055

Device Name: Lava™ Frame, Lava™ Frame Shade

Indications For Use: The Lava™ system is intended for CAD/CAM fabrication of all-ceramic dental restorations.

The system is used for the manufacturing of inlays, onlays, veneers, crowns and bridges.

Lava™ Frame and Lava™ Frame Shade are intended for the manufacturing of abutments. The titanium connection for the abutment must meet the following dimensions:

- Overall cementation surface > 30 mm<sup>2</sup>
- Height of the head of the titanium interface from the shoulder > 2.8 mm

The following systems fulfill the above described specifications:

- Co. Alltec Dental GmbH: Camlog Titanium-base for Ceramic-abutment – Abutment  $\phi \geq 4.3$  mm
- Co. Dentsply Friadent GmbH: Friadent Cera Base
- Co. Neoss GmbH: Neo Link Neoss Mono Abutment Titanium; Neo Link Neoss Multi Abutment Titanium; Neo Link Neoss Mono Aesthetic Abutment Titanium; Neo Link Neoss Multi Aesthetic Abutment Titanium; Matrix Abutment Hex – Regular Mono Titan; Matrix Abutment Hex – Regular Multi Titan; Matrix Abutment Hex – Narrow Mono Titan; Matrix Abutment Hex – Narrow Multi Titan; Matrix Abutment C-Lect – Regular Mono Titan; Matrix Abutment, C-Lect – Regular Multi Titan; Matrix Abutment C-Lect – Narrow Mono Titan; Matrix Abutment C-Lect – Narrow Multi Titan; Matrix Abutment ST – Mono Titan; Matrix Abutment ST – Mono Titan


Prescription Use    
 (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use    
 (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)

Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

Page 1 of \_\_\_\_\_

510(k) Number: K072055

92