



May 27, 2014

FDA CDRH DMC

JUN 01 2015

Received

Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center - WO66-G609  
10903 New Hampshire Avenue  
Silver Spring, Maryland 20993-0002

Re: 510(k) Submission

Att.: Document Mail Clerk

This is to notify you of the intention by 3Shape A/S to manufacture and market the below device:

**Type of 510(k) submission:** Traditional

**Proprietary Name:** 3Shape Abutment Designer™ Software

**Device Common Name:** Abutment Designer

**510(k) submitter:** 3Shape A/S, Holmens Kanal 7, DK-1060 Copenhagen K, Tel.: +45 7027 2620, Fax.: +45 7027 2621

**Contact Person:** Ms. Hanne Nielsen, Regulatory Affairs Manager, Tel.: +45 7027 2620

**Establishment Registration Number:** 3005940400

**Manufacturers Registration No.:** 10023901

**Classification Name:** Endosseous Dental Implant Abutment

**Regulation No.:** 872.3630

**Classification:** Class II

**Panel:** Dental

**Product Code:** NHA

**Submission copies:** 2 (eCopy and paper). The eCopy is an exact duplicate of the paper copy.

**Promotional Material:** Labeling examples and preliminary copies of the user documentation are enclosed.

**Substantial Equivalence:** The 3Shape Abutment Designer™ Software is similar in design and function to the Sirona Dental CAD/CAM System (K100152).

122

**Basis for the Submission:** New device

**Prior Submissions:** K133457 WD001  
CPT1200130 – Letter from the Agency  
I120810 – Pre-submission

This 510(k) submission is based on K133457. It takes into consideration the deficiencies raised by the Agency in email dated December 9th, 2014. This submission only relates to the part of the software which falls under product code NHA. To meet the Agency's requirements, a new software release was required and the version is therefore changed from DS2012-1 to DS2015-1. For a reply to the Additional Information request #2, please refer to the document "Deficiencies Previous Submission" which can be found in VOL\_001.

### **Design and Use of the Device**

The device is intended for prescription use (21 CFR 801 Subpart D).

The device is software only.

Best regards,  
3Shape A/S



Hanne Nielsen  
Regulatory Affairs Manager

Form Approved OMB No. 0910-0511 Expiration Date: April 30, 2016 See Instructions for OMB Statement

<p><b>DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION MEDICAL DEVICE USER FEE COVER SHEET</b></p>	<p><b>PAYMENT IDENTIFICATION NUMBER:</b> (b)(4) Write the Payment Identification number on your check.</p>
<p>A completed cover sheet must accompany each original application or supplement subject to fees. If payment is sent by U.S. mail or courier, please include a copy of this completed form with payment. Payment and mailing instructions can be found at: <a href="http://www.fda.gov/oc/mdufma/coversheet.html">http://www.fda.gov/oc/mdufma/coversheet.html</a></p>	
<p><b>1. COMPANY NAME AND ADDRESS (include name, street address, city state, country, and post office code)</b></p> <p>3SHAPE AS Holmens Kanal 7, floors 4</p> <p>KOBENHAVN 1060 DK</p> <p><b>1.1 EMPLOYER IDENTIFICATION NUMBER (EIN)</b></p>	<p><b>2. CONTACT NAME</b> Hanne Nielsen</p> <p><b>2.1 E-MAIL ADDRESS</b> hanne.nielsen@3shape.com</p> <p><b>2.2 TELEPHONE NUMBER (include Area code)</b> 45-70272620</p> <p><b>2.3 FACSIMILE (FAX) NUMBER (Include Area code)</b></p>
<p><b>3. TYPE OF PREMARKET APPLICATION (Select one of the following in each column; if you are unsure, please refer to the application descriptions at the following web site: <a href="http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm345263.htm">http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm345263.htm</a>)</b></p> <p><u>Select an application type:</u></p> <p><input checked="" type="checkbox"/> Premarket notification(510(k)); except for third party  <input type="checkbox"/> 513(g) Request for Information  <input type="checkbox"/> Biologics License Application (BLA)  <input type="checkbox"/> Premarket Approval Application (PMA)  <input type="checkbox"/> Modular PMA  <input type="checkbox"/> Product Development Protocol (PDP)  <input type="checkbox"/> Premarket Report (PMR)  <input type="checkbox"/> 30-Day Notice</p> <p><u>3.1 Select a center</u>  <input checked="" type="checkbox"/> CDRH  <input type="checkbox"/> CBER</p> <p><u>3.2 Select one of the types below</u>  <input checked="" type="checkbox"/> Original Application  <u>Supplement Types:</u>  <input type="checkbox"/> Efficacy (BLA)  <input type="checkbox"/> Panel Track (PMA, PMR, PDP)  <input type="checkbox"/> Real-Time (PMA, PMR, PDP)  <input type="checkbox"/> 180-day (PMA, PMR, PDP)</p>	
<p><b>4. ARE YOU A SMALL BUSINESS? (See the instructions for more information on determining this status)</b></p> <p><input type="checkbox"/> YES, I meet the small business criteria and have submitted the required qualifying documents to FDA  <input checked="" type="checkbox"/> NO, I am not a small business</p> <p><b>4.1 If Yes, please enter your Small Business Decision Number:</b></p>	
<p><b>5. FDA WILL NOT ACCEPT YOUR SUBMISSION IF YOUR COMPANY HAS NOT PAID AN ESTABLISHMENT REGISTRATION FEE THAT IS DUE TO FDA. HAS YOUR COMPANY PAID ALL ESTABLISHMENT REGISTRATION FEES THAT ARE DUE TO FDA?</b></p> <p><input checked="" type="checkbox"/> YES (All of our establishments have registered and paid the fee, or this is our first device, and we will register and pay the fee within 30 days of FDA's approval/clearance of this device.)  <input type="checkbox"/> NO (If "NO," FDA will not accept your submission until you have paid all fees due to FDA. This submission will not be processed; see <a href="http://www.fda.gov/cdrh/mdufma">http://www.fda.gov/cdrh/mdufma</a> for additional information)</p>	
<p><b>6. IS THIS PREMARKET APPLICATION COVERED BY ANY OF THE FOLLOWING USER FEE EXCEPTIONS? IF SO, CHECK THE APPLICABLE EXCEPTION.</b></p> <p><input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates  <input type="checkbox"/> This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only</p> <p><input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population  <input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially</p>	
<p><b>7. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION FOR WHICH FEES WERE WAIVED DUE TO SOLE USE IN A PEDIATRIC POPULATION THAT NOW PROPOSES CONDITION OF USE FOR ANY ADULT POPULATION? (If so, the application is subject to the fee that applies for an original premarket approval application (PMA).)</b></p>	

YES                       NO

**PAPERWORK REDUCTION ACT STATEMENT**

Public reporting burden for this collection of information is estimated to average 18 minutes per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to the address below.

Department of Health and Human Services, Food and Drug Administration, Office of Chief Information Officer, 8455 Colesville Road, COLE-14-14253 Silver Spring, MD 20993-0002

[Please do NOT return this form to the above address, except as it pertains to comments on the burden estimate.]

**8. USER FEE PAYMENT AMOUNT SUBMITTED FOR THIS PREMARKET APPLICATION**

(b)(4)

27-Feb-2015

Form FDA 3601 (05/13)

["Close Window"](#) [Print Cover sheet](#)

Form Approved OMB No. 0910-0511 Expiration Date: April 30, 2016 See Instructions for OMB Statement

<p><b>DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION MEDICAL DEVICE USER FEE COVER SHEET</b></p>	<p><b>PAYMENT IDENTIFICATION NUMBER:</b> (b)(4) Write the Payment Identification number on your check.</p>
<p>A completed cover sheet must accompany each original application or supplement subject to fees. If payment is sent by U.S. mail or courier, please include a copy of this completed form with payment. Payment and mailing instructions can be found at: <a href="http://www.fda.gov/oc/mdufma/coversheet.html">http://www.fda.gov/oc/mdufma/coversheet.html</a></p>	
<p><b>1. COMPANY NAME AND ADDRESS (include name, street address, city state, country, and post office code)</b></p> <p>3SHAPE AS Holmens Kanal 7, floors 4</p> <p>KOBENHAVN 1060 DK</p> <p><b>1.1 EMPLOYER IDENTIFICATION NUMBER (EIN)</b></p>	<p><b>2. CONTACT NAME</b> Hanne Nielsen</p> <p><b>2.1 E-MAIL ADDRESS</b> hanne.nielsen@3shape.com</p> <p><b>2.2 TELEPHONE NUMBER (include Area code)</b> 45-70272620</p> <p><b>2.3 FACSIMILE (FAX) NUMBER (Include Area code)</b></p>
<p><b>3. TYPE OF PREMARKET APPLICATION (Select one of the following in each column; if you are unsure, please refer to the application descriptions at the following web site: <a href="http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm345263.htm">http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm345263.htm</a>)</b></p> <p><u>Select an application type:</u></p> <p><input checked="" type="checkbox"/> Premarket notification(510(k)); except for third party  <input type="checkbox"/> 513(g) Request for Information  <input type="checkbox"/> Biologics License Application (BLA)  <input type="checkbox"/> Premarket Approval Application (PMA)  <input type="checkbox"/> Modular PMA  <input type="checkbox"/> Product Development Protocol (PDP)  <input type="checkbox"/> Premarket Report (PMR)  <input type="checkbox"/> 30-Day Notice</p> <p><u>3.1 Select a center</u>  <input checked="" type="checkbox"/> CDRH  <input type="checkbox"/> CBER</p> <p><u>3.2 Select one of the types below</u>  <input checked="" type="checkbox"/> Original Application  <u>Supplement Types:</u>  <input type="checkbox"/> Efficacy (BLA)  <input type="checkbox"/> Panel Track (PMA, PMR, PDP)  <input type="checkbox"/> Real-Time (PMA, PMR, PDP)  <input type="checkbox"/> 180-day (PMA, PMR, PDP)</p>	
<p><b>4. ARE YOU A SMALL BUSINESS? (See the instructions for more information on determining this status)</b></p> <p><input type="checkbox"/> YES, I meet the small business criteria and have submitted the required qualifying documents to FDA  <input checked="" type="checkbox"/> NO, I am not a small business</p> <p><b>4.1 If Yes, please enter your Small Business Decision Number:</b></p>	
<p><b>5. FDA WILL NOT ACCEPT YOUR SUBMISSION IF YOUR COMPANY HAS NOT PAID AN ESTABLISHMENT REGISTRATION FEE THAT IS DUE TO FDA. HAS YOUR COMPANY PAID ALL ESTABLISHMENT REGISTRATION FEES THAT ARE DUE TO FDA?</b></p> <p><input checked="" type="checkbox"/> YES (All of our establishments have registered and paid the fee, or this is our first device, and we will register and pay the fee within 30 days of FDA's approval/clearance of this device.)  <input type="checkbox"/> NO (If "NO," FDA will not accept your submission until you have paid all fees due to FDA. This submission will not be processed; see <a href="http://www.fda.gov/cdrh/mdufma">http://www.fda.gov/cdrh/mdufma</a> for additional information)</p>	
<p><b>6. IS THIS PREMARKET APPLICATION COVERED BY ANY OF THE FOLLOWING USER FEE EXCEPTIONS? IF SO, CHECK THE APPLICABLE EXCEPTION.</b></p> <p><input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates  <input type="checkbox"/> This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only</p> <p><input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population  <input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially</p>	
<p><b>7. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION FOR WHICH FEES WERE WAIVED DUE TO SOLE USE IN A PEDIATRIC POPULATION THAT NOW PROPOSES CONDITION OF USE FOR ANY ADULT POPULATION? (If so, the application is subject to the fee that applies for an original premarket approval application (PMA).)</b></p>	

YES  NO

**PAPERWORK REDUCTION ACT STATEMENT**

Public reporting burden for this collection of information is estimated to average 18 minutes per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to the address below.

Department of Health and Human Services, Food and Drug Administration, Office of Chief Information Officer, 8455 Colesville Road, COLE-14-14253 Silver Spring, MD 20993-0002

[Please do NOT return this form to the above address, except as it pertains to comments on the burden estimate.]

**8. USER FEE PAYMENT AMOUNT SUBMITTED FOR THIS PREMARKET APPLICATION**

(b)(4)

27-Feb-2015

Form FDA 3601 (05/13)

["Close Window"](#) [Print Cover sheet](#)

Form Approved OMB No. 0910-0511 Expiration Date: April 30, 2016 See Instructions for OMB Statement

<p><b>DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION MEDICAL DEVICE USER FEE COVER SHEET</b></p>	<p><b>PAYMENT IDENTIFICATION NUMBER:</b> (b)(4) Write the Payment Identification number on your check.</p>
<p>A completed cover sheet must accompany each original application or supplement subject to fees. If payment is sent by U.S. mail or courier, please include a copy of this completed form with payment. Payment and mailing instructions can be found at: <a href="http://www.fda.gov/oc/mdufma/coversheet.html">http://www.fda.gov/oc/mdufma/coversheet.html</a></p>	
<p><b>1. COMPANY NAME AND ADDRESS (include name, street address, city state, country, and post office code)</b></p> <p>3SHAPE AS Holmens Kanal 7, floors 4</p> <p>KOBENHAVN 1060 DK</p> <p><b>1.1 EMPLOYER IDENTIFICATION NUMBER (EIN)</b></p>	<p><b>2. CONTACT NAME</b> Hanne Nielsen</p> <p><b>2.1 E-MAIL ADDRESS</b> hanne.nielsen@3shape.com</p> <p><b>2.2 TELEPHONE NUMBER (include Area code)</b> 45-70272620</p> <p><b>2.3 FACSIMILE (FAX) NUMBER (Include Area code)</b></p>
<p><b>3. TYPE OF PREMARKET APPLICATION (Select one of the following in each column; if you are unsure, please refer to the application descriptions at the following web site: <a href="http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm345263.htm">http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm345263.htm</a>)</b></p> <p><u>Select an application type:</u></p> <p><input checked="" type="checkbox"/> Premarket notification(510(k)); except for third party  <input type="checkbox"/> 513(g) Request for Information  <input type="checkbox"/> Biologics License Application (BLA)  <input type="checkbox"/> Premarket Approval Application (PMA)  <input type="checkbox"/> Modular PMA  <input type="checkbox"/> Product Development Protocol (PDP)  <input type="checkbox"/> Premarket Report (PMR)  <input type="checkbox"/> 30-Day Notice</p> <p><u>3.1 Select a center</u>  <input checked="" type="checkbox"/> CDRH  <input type="checkbox"/> CBER</p> <p><u>3.2 Select one of the types below</u>  <input checked="" type="checkbox"/> Original Application  <u>Supplement Types:</u>  <input type="checkbox"/> Efficacy (BLA)  <input type="checkbox"/> Panel Track (PMA, PMR, PDP)  <input type="checkbox"/> Real-Time (PMA, PMR, PDP)  <input type="checkbox"/> 180-day (PMA, PMR, PDP)</p>	
<p><b>4. ARE YOU A SMALL BUSINESS? (See the instructions for more information on determining this status)</b></p> <p><input type="checkbox"/> YES, I meet the small business criteria and have submitted the required qualifying documents to FDA  <input checked="" type="checkbox"/> NO, I am not a small business</p> <p><b>4.1 If Yes, please enter your Small Business Decision Number:</b></p>	
<p><b>5. FDA WILL NOT ACCEPT YOUR SUBMISSION IF YOUR COMPANY HAS NOT PAID AN ESTABLISHMENT REGISTRATION FEE THAT IS DUE TO FDA. HAS YOUR COMPANY PAID ALL ESTABLISHMENT REGISTRATION FEES THAT ARE DUE TO FDA?</b></p> <p><input checked="" type="checkbox"/> YES (All of our establishments have registered and paid the fee, or this is our first device, and we will register and pay the fee within 30 days of FDA's approval/clearance of this device.)  <input type="checkbox"/> NO (If "NO," FDA will not accept your submission until you have paid all fees due to FDA. This submission will not be processed; see <a href="http://www.fda.gov/cdrh/mdufma">http://www.fda.gov/cdrh/mdufma</a> for additional information)</p>	
<p><b>6. IS THIS PREMARKET APPLICATION COVERED BY ANY OF THE FOLLOWING USER FEE EXCEPTIONS? IF SO, CHECK THE APPLICABLE EXCEPTION.</b></p> <p><input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates  <input type="checkbox"/> This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only</p> <p><input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population  <input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially</p>	
<p><b>7. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION FOR WHICH FEES WERE WAIVED DUE TO SOLE USE IN A PEDIATRIC POPULATION THAT NOW PROPOSES CONDITION OF USE FOR ANY ADULT POPULATION? (If so, the application is subject to the fee that applies for an original premarket approval application (PMA).)</b></p>	

YES  NO

**PAPERWORK REDUCTION ACT STATEMENT**

Public reporting burden for this collection of information is estimated to average 18 minutes per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to the address below.

Department of Health and Human Services, Food and Drug Administration, Office of Chief Information Officer, 8455 Colesville Road, COLE-14-14253 Silver Spring, MD 20993-0002

[Please do NOT return this form to the above address, except as it pertains to comments on the burden estimate.]

**8. USER FEE PAYMENT AMOUNT SUBMITTED FOR THIS PREMARKET APPLICATION**

(b)(4) 

27-Feb-2015

Form FDA 3601 (05/13)

["Close Window"](#) [Print Cover sheet](#)



Form Approved OMB No. 0910-0511 Expiration Date: April 30, 2016 See Instructions for OMB Statement

<p><b>DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION MEDICAL DEVICE USER FEE COVER SHEET</b></p>	<p><b>PAYMENT IDENTIFICATION NUMBER:</b> (b)(4) Write the Payment Identification number on your check.</p>																				
<p>A completed cover sheet must accompany each original application or supplement subject to fees. If payment is sent by U.S. mail or courier, please include a copy of this completed form with payment. Payment and mailing instructions can be found at: <a href="http://www.fda.gov/oc/mdufma/coversheet.html">http://www.fda.gov/oc/mdufma/coversheet.html</a></p>																					
<p><b>1. COMPANY NAME AND ADDRESS</b> (include name, street address, city state, country, and post office code)</p> <p>3SHAPE AS Holmens Kanal 7, floors 4</p> <p>KOBENHAVN 1060 DK</p> <p><b>1.1 EMPLOYER IDENTIFICATION NUMBER (EIN)</b></p>	<p><b>2. CONTACT NAME</b> Hanne Nielsen</p> <p><b>2.1 E-MAIL ADDRESS</b> hanne.nielsen@3shape.com</p> <p><b>2.2 TELEPHONE NUMBER</b> (include Area code) 45-70272620</p> <p><b>2.3 FACSIMILE (FAX) NUMBER</b> (Include Area code)</p>																				
<p><b>3. TYPE OF PREMARKET APPLICATION</b> (Select one of the following in each column; if you are unsure, please refer to the application descriptions at the following web site: <a href="http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm345263.htm">http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm345263.htm</a>)</p> <p>Select an application type:</p> <table border="0"> <tr> <td><input checked="" type="checkbox"/> Premarket notification(510(k)); except for third party</td> <td><b>3.1 Select a center</b></td> </tr> <tr> <td><input type="checkbox"/> 513(g) Request for Information</td> <td><input checked="" type="checkbox"/> CDRH</td> </tr> <tr> <td><input type="checkbox"/> Biologics License Application (BLA)</td> <td><input type="checkbox"/> CBER</td> </tr> <tr> <td><input type="checkbox"/> Premarket Approval Application (PMA)</td> <td><b>3.2 Select one of the types below</b></td> </tr> <tr> <td><input type="checkbox"/> Modular PMA</td> <td><input checked="" type="checkbox"/> Original Application</td> </tr> <tr> <td><input type="checkbox"/> Product Development Protocol (PDP)</td> <td><b>Supplement Types:</b></td> </tr> <tr> <td><input type="checkbox"/> Premarket Report (PMR)</td> <td><input type="checkbox"/> Efficacy (BLA)</td> </tr> <tr> <td><input type="checkbox"/> 30-Day Notice</td> <td><input type="checkbox"/> Panel Track (PMA, PMR, PDP)</td> </tr> <tr> <td></td> <td><input type="checkbox"/> Real-Time (PMA, PMR, PDP)</td> </tr> <tr> <td></td> <td><input type="checkbox"/> 180-day (PMA, PMR, PDP)</td> </tr> </table>		<input checked="" type="checkbox"/> Premarket notification(510(k)); except for third party	<b>3.1 Select a center</b>	<input type="checkbox"/> 513(g) Request for Information	<input checked="" type="checkbox"/> CDRH	<input type="checkbox"/> Biologics License Application (BLA)	<input type="checkbox"/> CBER	<input type="checkbox"/> Premarket Approval Application (PMA)	<b>3.2 Select one of the types below</b>	<input type="checkbox"/> Modular PMA	<input checked="" type="checkbox"/> Original Application	<input type="checkbox"/> Product Development Protocol (PDP)	<b>Supplement Types:</b>	<input type="checkbox"/> Premarket Report (PMR)	<input type="checkbox"/> Efficacy (BLA)	<input type="checkbox"/> 30-Day Notice	<input type="checkbox"/> Panel Track (PMA, PMR, PDP)		<input type="checkbox"/> Real-Time (PMA, PMR, PDP)		<input type="checkbox"/> 180-day (PMA, PMR, PDP)
<input checked="" type="checkbox"/> Premarket notification(510(k)); except for third party	<b>3.1 Select a center</b>																				
<input type="checkbox"/> 513(g) Request for Information	<input checked="" type="checkbox"/> CDRH																				
<input type="checkbox"/> Biologics License Application (BLA)	<input type="checkbox"/> CBER																				
<input type="checkbox"/> Premarket Approval Application (PMA)	<b>3.2 Select one of the types below</b>																				
<input type="checkbox"/> Modular PMA	<input checked="" type="checkbox"/> Original Application																				
<input type="checkbox"/> Product Development Protocol (PDP)	<b>Supplement Types:</b>																				
<input type="checkbox"/> Premarket Report (PMR)	<input type="checkbox"/> Efficacy (BLA)																				
<input type="checkbox"/> 30-Day Notice	<input type="checkbox"/> Panel Track (PMA, PMR, PDP)																				
	<input type="checkbox"/> Real-Time (PMA, PMR, PDP)																				
	<input type="checkbox"/> 180-day (PMA, PMR, PDP)																				
<p><b>4. ARE YOU A SMALL BUSINESS?</b> (See the instructions for more information on determining this status)</p> <p><input type="checkbox"/> YES, I meet the small business criteria and have submitted the required qualifying documents to FDA <input checked="" type="checkbox"/> NO, I am not a small business</p> <p><b>4.1</b> If Yes, please enter your Small Business Decision Number:</p>																					
<p><b>5. FDA WILL NOT ACCEPT YOUR SUBMISSION IF YOUR COMPANY HAS NOT PAID AN ESTABLISHMENT REGISTRATION FEE THAT IS DUE TO FDA. HAS YOUR COMPANY PAID ALL ESTABLISHMENT REGISTRATION FEES THAT ARE DUE TO FDA?</b></p> <p><input checked="" type="checkbox"/> YES (All of our establishments have registered and paid the fee, or this is our first device, and we will register and pay the fee within 30 days of FDA's approval/clearance of this device.)</p> <p><input type="checkbox"/> NO (If "NO," FDA will not accept your submission until you have paid all fees due to FDA. This submission will not be processed; see <a href="http://www.fda.gov/cdrh/mdufma">http://www.fda.gov/cdrh/mdufma</a> for additional information)</p>																					
<p><b>6. IS THIS PREMARKET APPLICATION COVERED BY ANY OF THE FOLLOWING USER FEE EXCEPTIONS? IF SO, CHECK THE APPLICABLE EXCEPTION.</b></p> <table border="0"> <tr> <td><input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates</td> <td><input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population</td> </tr> <tr> <td><input type="checkbox"/> This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only</td> <td><input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially</td> </tr> </table>		<input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates	<input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population	<input type="checkbox"/> This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only	<input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially																
<input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates	<input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population																				
<input type="checkbox"/> This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only	<input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially																				
<p><b>7. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION FOR WHICH FEES WERE WAIVED DUE TO SOLE USE IN A PEDIATRIC POPULATION THAT NOW PROPOSES CONDITION OF USE FOR ANY ADULT POPULATION?</b> (If so, the application is subject to the fee that applies for an original premarket approval application (PMA).)</p>																					

YES  NO

**PAPERWORK REDUCTION ACT STATEMENT**

Public reporting burden for this collection of information is estimated to average 18 minutes per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to the address below.

Department of Health and Human Services, Food and Drug Administration, Office of Chief Information Officer, 8455 Colesville Road, COLE-14-14253 Silver Spring, MD 20993-0002

[Please do NOT return this form to the above address, except as it pertains to comments on the burden estimate.]

**8. USER FEE PAYMENT AMOUNT SUBMITTED FOR THIS PREMARKET APPLICATION**

(b)(4) 

27-Feb-2015

Form FDA 3601 (05/13)

["Close Window"](#) [Print Cover sheet](#)

**Online Payment**

**Step 3: Confirm Payment**

1 | 2 | 3

**Thank you.**

**Your transaction has been successfully completed.**

(b)(4)



DEPARTMENT OF HEALTH AND HUMAN SERVICES  
Food and Drug Administration

Form Approved: OMB No. 0910-0120  
Expiration Date: January 31, 2017  
See PRA Statement below.

## Indications for Use

510(k) Number (if known)

Device Name

3Shape Abutment Designer™ Software

Indications for Use (Describe)

The 3Shape Abutment Designer™ Software is intended for designing endosseous dental implant abutments based on digital scans of 3D models.

Intended users are Dental Practitioners and Dental Laboratory staff.

Intended Operational Environment is Dental Laboratories.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

May 27, 2014

Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center - WO66-G609  
10903 New Hampshire Avenue  
Silver Spring, Maryland 20993-0002

Re: 510(k) Submission

Att.: Document Mail Clerk

This is to notify you of the intention by 3Shape A/S to manufacture and market the below device:

**Type of 510(k) submission:** Traditional

**Proprietary Name:** 3Shape Abutment Designer™ Software

**Device Common Name:** Abutment Designer

**510(k) submitter:** 3Shape A/S, Holmens Kanal 7, DK-1060 Copenhagen K, Tel.: +45 7027 2620, Fax.: +45 7027 2621

**Contact Person:** Ms. Hanne Nielsen, Regulatory Affairs Manager, Tel.: +45 7027 2620

**Establishment Registration Number:** 3005940400

**Manufacturers Registration No.:** 10023901

**Classification Name:** Endosseous Dental Implant Abutment

**Regulation No.:** 872.3630

**Classification:** Class II

**Panel:** Dental

**Product Code:** NHA

**Submission copies:** 2 (eCopy and paper). The eCopy is an exact duplicate of the paper copy.

**Promotional Material:** Labeling examples and preliminary copies of the user documentation are enclosed.

**Substantial Equivalence:** The 3Shape Abutment Designer™ Software is similar in design and function to the Sirona Dental CAD/CAM System (K100152).

**Basis for the Submission:** New device

**Prior Submissions:** K133457 WD001  
CPT1200130 – Letter from the Agency  
I120810 – Pre-submission


This 510(k) submission is based on K133457. It takes into consideration the deficiencies raised by the Agency in email dated December 9th, 2014. This submission only relates to the part of the software which falls under product code NHA. To meet the Agency's requirements, a new software release was required and the version is therefore changed from DS2012-1 to DS2015-1. For a reply to the Additional Information request #2, please refer to the document "Deficiencies Previous Submission" which can be found in VOL\_001.

**Design and Use of the Device**

The device is intended for prescription use (21 CFR 801 Subpart D).

The device is software only.

Best regards,  
3Shape A/S



Hanne Nielsen  
Regulatory Affairs Manager

## TABLE OF CONTENTS

### **ADMINISTRATIVE DOCUMENTS VOL 001**

<b>MEDICAL DEVICE USER FEE COVER SHEET AND INSTRUCTIONS</b>		
<b>PAYMENT AUTHORIZATION</b>		
<b>SUMBISSION COVER SHEET</b>		
<b>INDICATIONS FOR USE STATEMENT</b>		
<b>COVER LETTER</b>		<b>1-1</b>
<b>TABLE OF CONTENTS</b>		<b>1-3</b>
<b>510K SUMMARY</b>		<b>1-5</b>
<b>TRUTHFUL AND ACCURATE STATEMENT</b>		<b>1-9</b>
<b>PROPOSED LABELLING</b>		<b>1-10</b>
<b>SPECIFICATIONS</b>		<b>1-11</b>
<b>SUBSTANTIAL EQUIVALENCE COMPARISON</b>		<b>1-12</b>
<b>BIOCOMPATIBILITY</b>		<b>1-15</b>
<b>EMC AND ELECTRICAL SAFETY</b>		<b>1-16</b>
<b>SHELF LIFE</b>		<b>1-17</b>
<b>STERILIZATION</b>		<b>1-18</b>
<b>CORROSION TESTING</b>		<b>1-19</b>
<b>IMPLANT TO ABUTMENT COMPATIBILITY</b>		<b>1-20</b>
<b>MATERIAL COMPOSITION</b>		<b>1-21</b>
<b>MECHANICAL PROPERTIES</b>		<b>1-22</b>
<b>MODIFIED SURFACES INFORMATION</b>		<b>1-23</b>
<b>DEFICIENCIES PREVIOUS SUBMISSION</b>		<b>1-24</b>
<b>APPLICABLE VOLUNTARY STANDARDS</b>		<b>1-27</b>
510K ACCEPTANCE CHECKLIST	APPENDIX 1 1 (26 PAGES)	
STANDARD REPORTS – FORM 3654	APPENDIX 1 2 (9 PAGES)	

### **LEVEL OF CONCERN VOL 002**

<b>LEVEL OF CONCERN</b>	<b>2-1</b>
-------------------------	------------

### **SOFTWARE DESCRIPTION VOL 003**

<b>SOFTWARE DESCRIPTION</b>	<b>3-1</b>
-----------------------------	------------

### **DEVICE HAZARD ANALYSIS VOL 004**

<b>DEVICE HAZARD ANALYSIS</b>	<b>4-1</b>
ABUTMENT DESIGNER PRODUCT RISK ASSESSMENT	APPENDIX 4 1 (1 DOCUMENT)

### **SOFTWARE REQUIREMENTS SPECIFICATION VOL 005**

<b>SOFTWARE REQUIREMENTS SPECIFICATION</b>	<b>5-1</b>
REQUIREMENTS SPECIFICATION DS2015 1 RSXXXX	APPENDIX 5 1 (8 DOCUMENTS)



**ARCHITECTURE DESIGN CHART** **VOL 006**

**ARCHITECTURE DESIGN CHART** **6-1**

**SOFTWARE DESIGN SPECIFICATION** **VOL 007**

**SOFTWARE DESIGN SPECIFICATION** **7-1**  
 SOFTWARE SPECIFICATION DS2015 1 SSXXXX APPENDIX 7 1 (9 DOCUMENTS)

**TRACEABILITY ANALYSIS** **VOL 008**

**TRACEABILITY ANALYSIS** **8-1**

**SOFTWARE DEVELOPMENT ENVIRONMENT DESCRIPTION** **VOL 009**

**SOFTWARE DEVELOPMENT ENVIRONMENT DESCRIPTION** **9-1**

**VERIFICATION AND VALIDATION DOCUMENTATION** **VOL 010**

**VERIFICATION AND VALIDATION DOCUMENTATION** **10-1**  
**ABUTMENT DESIGNER BETA REPORT – ABSTRACT** **10-3**  
 RELEASE PROTOCOL D2015-1-RP-2.15.2.0 APPENDIX 10 1 (1 DOCUMENT)  
 TEST RESULTS DS2015 1 RSXXXX TR APPENDIX 10 2 (8 DOCUMENTS)  
 TEST RESULTS DS2015 1 SSXXXX TR APPENDIX 10 3 (10 DOCUMENTS)

**REVISION LEVEL HISTORY** **VOL 011**

**REVISION LEVEL HISTORY** **11-1**

**UNRESOLVED ANOMALIES** **VOL 012**

**UNRESOLVED ANOMALIES** **12-1**

**LABELING** **VOL 013**

**LABELING** **13-1**  
 3SHAPE DENTAL SYSTEM SOFTWARE BROCHURE (EXCERPT) APPENDIX 13 1 (5 PAGES)  
 MANUAL DS-2.15.2.0-A-EN VERSION (EXCERPT) APPENDIX 13 2 (62 PAGES)



3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

**510(K) SUMMARY – Traditional 510(K)**

**Submitter Information**

A Company Name: 3Shape A/S  
B Company Address: Holmens Kanal 7  
DK-1060 Copenhagen K  
C Company Phone: +45 7027 2620  
Company Fax: +45 7027 2621  
D Contact Person: Hanne Nielsen  
Regulatory Affairs Manager  
E Date Summary Prepared: May 27, 2015

**Device Identification**

A Trade/proprietary Name: 3Shape Abutment Designer™  
Software  
B Common Name: Abutment Designer  
C Device Classification Name: Endosseous Dental Implant Abutment  
C Regulation Number: 872.3630  
C Classification: Class II  
D Product Code: NHA

**Predicate Device**

The 3Shape Abutment Designer™ Software has equivalent intended use and technical characteristics as the Sirona Dental CAD/CAM System (K100152).

3Shape Abutment Designer™ Software does not provide any physical parts that can come into contact with the patient.

Therefore, the differences between the Device and the predicates do not raise additional concerns concerning the Device's safety and effectiveness.

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

### **Intended Use**

The 3Shape Abutment Designer™ Software is intended for designing endosseous dental implant abutments based on digital scans of 3D models.

### **Device Description**

The 3Shape Abutment Designer™ Software receives scan data containing topographical characteristics of real teeth, dental impressions, or stone models. The software provides the user with the ability to create matching endosseous dental implant abutments using computer aided design.

The output of the device is a computer file containing the abutment(s) in digital form which can be used by 3rd party providers of milling machines and/or rapid prototyping machines to manufacture physical abutments.

The 3Shape Abutment Designer™ Software can be run on properly configured "off-the-shelf" PC hardware running Microsoft® Windows and using a standard keyboard and mouse.

### **Scientific Concept**

The underlying scientific concept is the use of CAD/CAM technology for the optical acquisition of the topographical characteristics of dental impressions, and models and the design of individual mesostructures using recorded data (CAD).

### **Materials Used**

Software. Not applicable.

### **Physical Properties**

Software. Not applicable.

## 3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

**Summary of the technological characteristics**

The 3Shape Abutment Designer™ is a software only device programmed in (b)(4) and has the following PC/laptop requirements:

<b>Item</b>	<b>Minimum Requirements</b>	<b>Recommended</b>
<b>OS</b>	Windows 7 32-bit Professional*	Windows 7 64-bit Professional
<b>RAM</b>	4GB	8GB (16GB)
<b>Video Card</b>	512MB DirectX 10 (1GB DirectX 10) NVIDIA GeForce	1GB DirectX 11 (2GB DirectX 11) NVIDIA GeForce
<b>Available HDD Space</b>	250GB	500GB (1TB)
<b>CPU</b>	Intel Core i5 or equivalent	Intel Core i7 or equivalent
<b>Monitor resolution</b>	1440 x 900 pixels	1920 x 1080 pixels
<b>3D Mouse</b>	None	3DConnexion SpaceMouse™ Pro
<b>Network</b>	Internet connection	
<b>USB ports</b>	USB 2.0 port for 3Shape desktop scanner	
<b>Mouse</b>	Mouse with wheel button support	

**Nonclinical Testing**

Software, hardware, and integration verification and validation testing was performed in accordance with the FDA Guidance Document "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices" (Issued on May 11, 2005).

Prior to release, verification and validation testing of the 3Shape Abutment Designer™ Software has been completed using the approved acceptance criteria: Each user need has its own validation acceptance criteria; each specification has its own verification acceptance criteria; bug verification consists in ensuring issue is not reproducible; issues reported by beta partners must be reviewed and handled appropriately; beta partners must report functionality to be an improvement over previous version.

The validation suite includes validation of implemented mitigations related to device hazards identified in the risk management procedures.

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

All test results have been reviewed and approved, showing the 3Shape Abutment Designer™ Software to be safe and effective.

**Clinical Testing**

Clinical testing is not a requirement and has not been performed.

**Conclusion**

Based on a comparison of intended use, indications, construction materials, principle of operations, features and technical data, and the test results, the 3Shape Abutment Designer™ Software is found to be as safe and effective as the Predicate Devices. Intended use and performance is found to be substantially equivalent to the Predicate Devices.



## Premarket Notification Truthful And Accurate Statement

[As Required by 21 CFR 807.87(k)]

I certify that, in my capacity as Regulatory Affairs Manager of 3Shape A/S, I believe to the best of my knowledge, that all data and information submitted in the premarket notification are truthful and accurate and that no material fact has been omitted.

A handwritten signature in blue ink, appearing to read "Hanne Nielsen", written over a horizontal line.

(Signature)

Hanne Nielsen

A handwritten date in blue ink, "MAY 27, 2015", written over a horizontal line.

(Date)

\_\_\_\_\_  
\*(Premarket Notification [510(k)] Number)

## 3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

**Proposed Labelling**

The 3Shape Abutment Designer™ Software is a prescription device and is exempt from needing adequate directions for lay use.

The following labelling is included (please see "VOL 013 Labelling" of this submission):

- PDF copy of User Manual
- Dental System™ Brochure

The application is for download only, and hence there is no CD label.

**Device Specific Requirements**

There are no applicable requirements in a device specific regulation.

**Special Controls Document**

The 3Shape Abutment Designer™ Software falls into the following product group:

**NHA**

The guidance document entitled "Class II Special Controls Guidance Document: Root Form Endosseous Dental Implants and Endosseous Dental Implant Abutments" serve as the special control.

The Guidance Document lists the below labelling requirements.

<b>Requirement</b>	<b>Referenced in</b>
Labelling in sufficient detail to satisfy the requirements of 21 CFR 807.87(e)	Online help Safety and Setup Guide
Provide users with a surgical manual along with the instructions for use	N/A The 3Shape Abutment Designer™ Software is not used in surgery
Provide all relevant precautions and warnings in the professional labelling	Online help Safety and Setup Guide
Precautions or warnings that relate to unpackaging or sterility	N/A The 3Shape Abutment Designer™ Software is a software device and not supplied sterile.
If any parts are provided non sterile we recommend that you provide sterilization instructions	N/A The 3Shape Abutment Designer™ Software is a software device and not supplied sterile.
If patient labelling is appropriate, we recommend that you follow Guidance on Medical Device Patient Labelling; Final Guidance for Industry and FDA Reviewers	N/A 3Shape Abutment Designer™ Software is a prescription device and patient labeling is not required.

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

---

## Specifications

### 1. Device Description

The 3Shape Abutment Designer™ Software falls under the device specific FDA guidance document "[Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices, May 11 2005](#)".

Thus, a narrative description of the Device, indications for use, principles of operation, power source, composition, and other information necessary to understand the device can be found under the Volume "VOL 003 Software Description" in the document "001 Software Description" included in this submission.

The 3Shape Abutment Designer™ Software does not include any physical part(s), therefore a physical description is not included in this submission.

## Substantial Equivalence Comparison

### 1. Predicates

The 3Shape Abutment Designer™ Software has the same intended uses and technical characteristics as the Sirona Dental CAD/CAM System (K100152) as listed in "Table 1: Predicate"

**Table 1: Predicate**

<b>Predicate</b>	<b>Manufacturer</b>	<b>510(k) number</b>	<b>Product code</b>
Dental CAD/CAM System	Sirona	K100152	NHA*

\* Endosseous dental implant abutments, 21CFR872.3630



## 2. Intended Use Comparison

### 2.1. 3Shape Abutment Designer™ Software

The Device's Intended Use, Intended users, and Intended Operational is stated in the "VOL 001 Administrative Documents" volume of this submission and reproduced here:

*"The 3Shape Abutment Designer™ Software is intended for designing endosseous dental implant abutments based on digital scans of 3D models.*

*Intended users are Dental Practitioners and Dental Laboratory staff.*

*Intended Operational Environment is Dental Laboratories."*

In the following sections, the similarities with the predicate are discussed.

### 2.2. Dental CAD/CAM System

The predicate's Intended Use can be extracted from the FDA 510(k) Premarket Notification Database:

"The Sirona Dental CAD/CAM System is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations. The system consists of three major parts: TiBase, InCoris mesostructure, and CAD/CAM software. Specifically, the InCoris mesostructure and TiBase components make up a two-piece abutment which is used in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity.

The InCoris mesostructure may also be used in conjunction with the Camlog Titanium base CAD/CAM (types K2244.xxxx) (K083496) in the Canmlog Implant System. The CAD/CAM software is intended to design and fabricate the InCoris mesostructure. The InCoris mesostructure and TiBase two-piece abutment is compatible with the following implants systems:

[List of compatible implant systems omitted for brevity]

Underlined segments indicate similarity with 3Shape Abutment Designer™ Software.

Note, the predicate is a CAD/CAM System bundled with physical two piece abutments; the 3Shape Abutment Designer™ Software is CAD/CAM only.

### 3. Characteristics Comparison

A Characteristics Comparison can be seen in "Table 2: Substantial Equivalence Chart".

**Table 2: Substantial Equivalence Chart**

<b>Feature name</b>	<b>3Shape Abutment Designer™ Software</b>	<b>Sirona Dental CAD/CAM System (K100152)</b>
Graphical UI	Yes	Yes
Windows OS platform	Yes	Yes
Uses standard PC hardware	Yes	Yes
Digitally imports topography of teeth by 3D Scan	Yes	Yes
Uses 3D CAD design tools	Yes	Yes
Custom abutment design	Yes	Yes
Screw retained design	Yes	Yes
Implant Bar design	Yes	Yes
Export to remote milling machine by internet	Yes	Yes
Network Protocol	Internet/TCP IP	Internet/TCP IP
Intended users	Dental practitioners and dental labs	Dental practitioners and dental labs
Output type	Computer file	Computer file
Device submission includes pre manufactured prosthetics*	No	Yes

\* Endosseous dental implant abutments as per 21CFR872.3630

### 4. Conclusion

The 3Shape Abutment Designer™ Software and the predicate only deviate significantly in the cases where the predicate is bundled with a physical dental implant abutment. 3Shape Abutment Designer™ Software does not provide any physical parts that can come into contact with the patient.

The differences between the 3Shape Abutment Designer™ Software and the predicate do not raise additional concerns with respect to the safety and effectiveness of the 3Shape Abutment Designer™ Software.

Based on the information presented, we conclude Substantial Equivalence between the predicate and the 3Shape Abutment Designer™ Software.

---

## **Biocompatibility**

### Biocompatibility

The 3Shape Abutment Designer™ Software has no patient contacting components and thus biocompatibility requirements do not apply.

---

## **EMC and Electrical Safety**

### EMC and Electrical Safety

The 3Shape Abutment Designer™ Software has no patient or user contacting components as it is a software device. EMC and Electrical Safety evaluation is therefore not applicable.

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

---

## Shelf Life

### Shelf Life

The 3Shape Abutment Designer™ Software is a software device and thus shelf life requirements do not apply.

---

## **Sterilization**

### Sterilization

The 3Shape Abutment Designer™ Software is a software device and thus sterilization requirements do not apply.

---

## **Corrosion Testing**

The 3Shape Abutment Designer™ Software is a software device only. Corrosion Testing is therefore not applicable.

---

## **Implant to Abutment Compatibility**

The 3Shape Abutment Designer™ Software is a software device only.

The output of the device is a computer file containing the abutment(s) in digital form which can be used by 3<sup>rd</sup> party providers of milling machines and/or rapid prototyping machines to manufacture physical abutments. The 3Shape Abutment Designer™ Software does not provide any means of creating physical abutments, nor does it consist of any physical parts.

A digital representation of the Implant is required for the software to work, but only to ensure the Abutment will fit on the Seating Geometry / Interface of the Implant. The 3Shape Abutment Designer™ does not provide any means to design, alter, or manufacture any part of the Implant including, but not limited to, the abutment to implant interface.

Correct compatibility with the digital representation of the Implant is validated as described in the "VOL 10 Verification and Validation" Volume of this submission.



3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

---

## Material Composition

### Material Identity

The 3Shape Abutment Designer™ Software is a software device only. Material Identity is therefore not applicable.

### Chemical Composition and anticipated impurities

The 3Shape Abutment Designer™ Software is a software device only. Chemical composition and anticipated impurities are therefore not applicable.

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

---

## **Mechanical Properties**

The 3Shape Abutment Designer™ Software is a software device only. Mechanical Properties are therefore not applicable.

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

---

## **Modified Surfaces Information**

The 3Shape Abutment Designer™ Software is a software device only. Modified Surfaces Information is therefore not applicable.

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

---

## Deficiencies Previous Submission

(b)(4)



3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

(b)(4)



3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

(b)(4)



---

## **Applicable Voluntary Standards**

### **A. Radiation control for Health and Safety Act**

- i. As the 3Shape Abutment Designer™ is a software only application this is not applicable.

### **B. CDRH – Recognized Voluntary Standards**

- i. There are not recognized voluntary standards applicable for the 3Shape Abutment Designer™

### **C. Standard Data Reports**

- i. A standard data report has been filled out for standards ISO 13485, ISO 14971 and IEC 62304, included as appendixes.
- ii. Summary report tables are included in the standard data reports for standards ISO 13485, ISO 14971 and IEC 62304 are included as appendixes.

Appendix - Acceptable Checklists

*Contains Nonbinding Recommendations*

**Acceptance Checklist  
for Traditional 510(k)s**

**(should be completed within 15 days of DCC receipt)**

*The following information is not intended to serve as a comprehensive review.*

**510(k) Number:** \_\_\_\_\_ **Date Received by DCC:** \_\_\_\_\_

**Lead Reviewer Name:** \_\_\_\_\_ **Branch:** \_\_\_\_\_ **Division:** \_\_\_\_\_ **Office:** \_\_\_\_\_

Note: If an element is left blank on the checklist, it does not mean the checklist is incomplete; it means the reviewer did not assess the element during RTA and that element will be assessed during substantive review.

Preliminary Questions		
Answers in the shaded blocks indicate consultation with Center advisor is needed.	Yes	No
<p><b>1. Is the product a device (per section 201(h) of the FD&amp;C Act) or a combination product (per 21 CFR 3.2(e)) with a device constituent part subject to review in a 510(k)?</b></p> <p>If it appears not to be a device (per section 201(h) of the FD&amp;C Act) or such a combination product, or you are unsure, consult with the CDRH Jurisdictional Officer or the CBER Office Jurisdiction Liaison to determine the appropriate action, and inform division management. <i>Provide a summary of the Jurisdictional Officer's/Liaison's determination.</i> If the product does not appear to be a device or such a combination product, mark "No."</p>	X	
<p><b>Comments:</b></p>		
<p><b>2. Is the application with the appropriate Center?</b></p> <p>If the product is a device or a combination product with a device constituent part, is it subject to review by the Center in which the submission was received? If you believe the application is not with the appropriate Center or you are unsure, consult with the CDRH Jurisdictional Officer or CBER Office Jurisdiction Liaison to determine the appropriate action and inform your division management. <i>Provide a summary of the Jurisdictional Officer's/Liaison's determination.</i> If application should not be reviewed by your Center mark "No."</p>	X	
<p><b>Comments:</b></p>		
<p><b>3. If a Request for Designation (RFD) was submitted for the device or combination product with a device constituent part and assigned to your center, identify the RFD # and confirm the following:</b></p>	N/A	

Acceptance Checklist for Traditional 510(k)



*Contains Nonbinding Recommendations*

<p>a) <b>Is the device or combination product the same (e.g., design, formulation) as that presented in the RFD submission?</b></p> <p>b) <b>Are the indications for use for the device or combination product identified in the 510(k) the same as those identified in the RFD submission?</b></p> <p>If you believe the product or the indications presented in the 510(k) have changed from the RFD, or you are unsure, consult with the CDRH Jurisdictional Officer or appropriate CBER Jurisdiction Liaison to determine the appropriate action and inform your division management. <i>Provide summary of Jurisdictional Officer's/Liaison's determination.</i></p> <p>If the answer to either question above is no, mark "No." If there was no RFD, skip this question.</p>		
<p><b>Comments:</b></p>		
<p><b>4. Is this device type eligible for a 510(k) submission?</b></p> <p>If a 510(k) does not appear to be appropriate (e.g., Class III type and PMA required, or Class I or II type and 510(k)-exempt), you should consult with the CDRH 510(k) Program Director or appropriate CBER staff during the acceptance review. If 510(k) is not the appropriate regulatory submission, mark "No."</p>	X	
<p><b>Comments:</b></p>		
<p><b>5. Is there a pending PMA for the same device with the same indications for use?</b></p> <p>If yes, consult division management and the CDRH 510(k) Program Director or appropriate CBER staff to determine the appropriate action.</p>		X
<p><b>Comments:</b></p>		
<p><b>6. If clinical studies have been submitted, is the submitter the subject of an Application Integrity Policy (AIP)?</b></p> <p>If yes, consult with the CDRH Office of Compliance/Division of Bioresearch Monitoring (OC/DBM - BIMO) or CBER Office of Compliance and Biologics Quality/Division of Inspections and Surveillance/Bioresearch Monitoring Branch (OCBQ/DIS/BMB) to determine the appropriate action. Check on web at <a href="http://www.fda.gov/ICECI/EnforcementActions/ApplicationIntegrityPolicy/ucm134453.htm">http://www.fda.gov/ICECI/EnforcementActions/ApplicationIntegrityPolicy/ucm134453.htm</a>.</p>		N/A

If the answer to 1 or 2 appears to be "No," then stop review of the 510(k) and issue the "Original Jurisdictional Product" letter.  
 If the answer to 3a or 3b appears to be "No," then stop the review and contact the CDRH Jurisdictional Officer or CBER Office of Jurisdiction Liaison.  
 If the answer to 4 is "No", the lead reviewer should consult division management and other Center resources to determine the appropriate action.  
 If the answer to 5 is "Yes," then stop review of the 510(k), contact the CDRH 510(k) Staff and PMA Staff, or appropriate CBER staff.

Acceptance Checklist for Traditional 510(k)

***Contains Nonbinding Recommendations***

If the answer to 6 is “Yes,” then contact CDRH/OC/DBM – BIMO or CBER/OCBQ/DIS/BMB, provide a summary of the discussion with the BIMO Staff, and indicate BIMO’s recommendation/action.

<b><u>Organizational Elements</u></b>		
<i>Failure to include these items alone generally should not result in an RTA designation</i>		
	<b>Yes</b>	<b>No</b>
a. Submission contains Table of Contents	<input checked="" type="checkbox"/>	<input type="checkbox"/>
b. Each section is labeled (e.g., headings or tabs designating Device Description section, Labeling section, etc.)	<input checked="" type="checkbox"/>	<input type="checkbox"/>
c. All pages of the submission are numbered <i>All pages should be numbered in such a manner that information can be referenced by page number. This may be done either by consecutively numbering the entire submission, or numbering the pages within a section (e.g., 12-1, 12-2...).</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
d. Type of 510(k) is identified– traditional, abbreviated, or special <i>If type of 510(k) is not designated, review as a traditional</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Comments:		

<b><u>Elements of a Complete Submission (RTA Items)</u></b>				
<b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>				
Submission should be designated RTA if not addressed				
<b>Check “Yes” if item is present, “N/A” if it is not needed and “No” if it is not included but needed.</b>				
	<ul style="list-style-type: none"> <li>Any “No” answer will result in a “Refuse to Accept” decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>	<b>Yes</b>	<b>N/A</b>	<b>No</b>
<b>A.</b>	<b>Administrative</b>			
	1. All content used to support the submission is written in English (including translations of test reports, literature articles, etc.)	<input checked="" type="checkbox"/>		<input type="checkbox"/>
	Comments:			
	2. Submission identifies the following (such as in CDRH Premarket Review Submission Cover Sheet ( <a href="#">Form 3514</a> ) or 510(k) cover letter):	<input checked="" type="checkbox"/>		<input type="checkbox"/>
	a. Device trade name or proprietary name	<input checked="" type="checkbox"/>		<input type="checkbox"/>
	b. Device common name	<input checked="" type="checkbox"/>		<input type="checkbox"/>

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>					
Submission should be designated RTA if not addressed					
<b>Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.</b>					
		<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>	<b>Yes</b>	<b>N/A</b>	<b>No</b>
	c.	Device class and panel or Classification regulation or Statement that device has not been classified with rationale for that conclusion	<input checked="" type="checkbox"/>		<input type="checkbox"/>
		Comments:			
	3.	Submission contains Indications for Use Statement with Rx and/or OTC designated (see also 21 CFR 801.109) <i>Submitter should use format appropriate for the reviewing Center/Office (CDRH/ODE, CDRH/OIVD, CBER/OBRR, CBER/OCTGT). If not provided in correct format, request the correct format during substantive review.</i>	<input checked="" type="checkbox"/>		<input type="checkbox"/>
		Comments:			
	4.	Submission contains 510(k) Summary or 510(k) Statement <i>Either a) or b) must be answered "Yes" to be considered complete. Identify any missing element(s) in Comments.</i>	<input checked="" type="checkbox"/>		<input type="checkbox"/>
	a.	Summary contains all elements per 21 CFR 807.92 <i>See also <a href="#">510(k) Summary Checklist</a></i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	b.	Statement contains all elements per 21 CFR 807.93	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Comments: VOL_001			
	5.	Submission contains Truthful and Accuracy Statement per 21 CFR 807.87(k) <i>See recommended <a href="#">format</a>. Select "Yes" if statement is present and includes the text in the recommended format, and is signed by a responsible person of the firm (not consultant).</i>	<input checked="" type="checkbox"/>		<input type="checkbox"/>
		Comments: VOL_001			

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>					
Submission should be designated RTA if not addressed					
Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.					
			Yes	N/A	No
		<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>			
6.		Submission contains Class III Summary and Certification <i>See recommended <a href="#">content</a>. Form should be signed by a responsible person of the firm, not a consultant. Select "N/A" only if submission is not a Class III 510(k).</i>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
		Comments:			
7.		Submission contains clinical data <i>Select "N/A" if the submission does not contain clinical data. If "N/A" is selected, parts a and b below are omitted from the checklist.</i>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
	a.	Submission includes completed Financial Certification ( <a href="#">FDA Form 3454</a> ) or Disclosure ( <a href="#">FDA Form 3455</a> ) information for each covered clinical study included in the submission. <i>Select "N/A" if the submitted clinical data is not a "covered clinical study" as defined in the <a href="#">Guidance for Industry- Financial Disclosures by Clinical Investigators</a></i>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	b.	Submission includes completed Certification of Compliance with requirements of ClinicalTrials.gov Data Bank ( <a href="#">FDA Form 3674</a> ) (42 U.S.C. 282(j)(5)(B)) for each applicable device clinical trial included in the submission. <i>Select "N/A" if the submitted clinical data is not an "applicable device clinical trial" as defined in <a href="#">Title VIII of FDAAA, Sec. 801(j)</a></i>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Comments:			
8.		If submission references use of a national or international standard as part of demonstration of substantial equivalence, submission contains complete Standards Data Report for 510(k)s ( <a href="#">FDA Form 3654</a> ) <i>There should be a completed form for each referenced national or international standard. Select "N/A" only if submission does not reference any standards.</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>				
Submission should be designated RTA if not addressed				
<b>Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.</b>				
	<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>	Yes	N/A	No
	Comments:			
9.	<p>The submission identifies prior submissions for the same device for which FDA provided feedback related to the data or information needed to support substantial equivalence (e.g., submission numbers for Pre-Submission, IDE, prior not substantially equivalent (NSE) determination, prior 510(k) that was deleted or withdrawn) or states that there were no prior submissions for the subject device.</p> <p><i>This information may be included in the Cover Letter (i.e., as a statement that there were no prior submissions for the device or a listing of the number(s) of the prior submissions). Alternatively, a list of submission numbers may be found in Section F (prior related submissions section) of the CDRH Coversheet form (Form 3514) to address this criterion. Please be advised that if this section of the form is left blank, it should not be considered a statement that there were no prior submissions.</i></p>	<input checked="" type="checkbox"/>		<input type="checkbox"/>
a.	<p>If there were prior submissions, the submitter has identified where in the current submission any issues related to a determination of substantial equivalence outlined in prior communications are addressed.</p> <p><i>To address this criterion, the submission may include a separate section with the prior submission number(s), a copy of the FDA feedback (e.g., letter, meeting minutes), and a statement of how or where in the submission this prior feedback was addressed. Note that the adequacy of how the feedback was addressed should be assessed during the substantive review. For additional information regarding the Pre-Submission process, please refer to the Draft Guidance "<a href="http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm310375.htm">Medical Devices: The Pre-Submission Program and Meetings with FDA Staff</a>." (<a href="http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm310375.htm">http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm310375.htm</a>). Once finalized, this guidance will represent the</i></p>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>									
Submission should be designated RTA if not addressed									
<b>Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.</b>									
	<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>				<b>Yes</b>	<b>N/A</b>	<b>No</b>		
		Agency's current thinking on this topic. <i>Select "N/A" if the submitter states there were no prior submissions in criterion above.</i>							
		Comments: DEFICIENCIES PREVIOUS SUBMISSION IN VOL - 001							
<b>B.</b>	<b>Device Description</b>								
	10.	a.	<p>If there are requirements regarding the device description, such as special controls, in a device-specific regulation that are applicable to the device, the submission includes device description information to establish that the submitter has followed the device-specific requirement.</p> <p><i>Select "N/A" if there are no applicable requirements in a device-specific regulation. Select "No" if the submission does not include a rationale for any omitted information. Note that the adequacy of how such requirements have been addressed should be assessed during the substantive review.</i></p>			<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		b.	<p>If there is a device-specific guidance, other than a special controls guidance document, applicable to the device, the submission includes device description information to establish that the submitter has addressed the recommendations or otherwise has met the applicable statutory or regulatory criteria through an alternative approach.</p> <p><i>Select "N/A" if there is no applicable device-specific guidance. Select "No" if the submission does not include a rationale for any omitted information or any alternative approach as outlined above. Note that the adequacy of how recommendations in a device-specific guidance, etc., have been addressed should be assessed during the substantive review.</i></p>			<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Comments:							

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>				
Submission should be designated RTA if not addressed				
Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.				
		Yes	N/A	No
		<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>		
11.	Descriptive information is present and consistent within the submission (e.g., the device description section is consistent with the device description in the labeling), including:			
a.	A description of the principle of operation and mechanism of action for achieving the intended effect.		<input checked="" type="checkbox"/>	<input type="checkbox"/>
b.	A description of proposed conditions of use, such as surgical technique for implants; anatomical location of use; user interface; how the device interacts with other devices; and/or how the device interacts with the patient.		<input checked="" type="checkbox"/>	<input type="checkbox"/>
c.	A list and description of each device for which clearance is requested. <i>Select "N/A" if there is only one device or model. "Device" may refer to models, part numbers, or various sizes, etc.</i>		<input type="checkbox"/>	<input checked="" type="checkbox"/>
Comments:				
12.	Submission contains representative engineering drawing(s), schematics, illustrations and/or figures of the device that are clear, legible, labeled, and include dimensions. <i>In lieu of drawings, schematics, etc. of each device to be marketed, "representative" drawings, etc. may be provided, where "representative" is intended to mean that the drawings, etc. provided capture the differences in design, size, and other important characteristics of the various models, sizes, or versions of the device(s) to be marketed.</i> <i>Select "N/A" if the submitter provided a rationale for why the submission does not contain engineering drawings, schematics, etc. (e.g., device is a reagent and figures are not pertinent to describe the device).</i>			<input type="checkbox"/>
Comments:				

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>					
Submission should be designated RTA if not addressed					
Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.					
			Yes	N/A	No
		<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>			
	13.	If device is intended to be marketed with multiple components, accessories, and/or as part of a system, <i>Select "N/A" if the device is not intended to be marketed with multiple components, accessories, and/or as part of a system.</i>		<input type="checkbox"/>	
		a. Submission includes a list of all components and accessories to be marketed with the subject device.	<input checked="" type="checkbox"/>		<input type="checkbox"/>
		b. Submission includes a description (as detailed in item 11.a. and b. and 12 above) of each component or accessory. <i>Select "N/A" if the component(s)/accessory(ies) has been previously cleared, or is exempt, and the proposed indications for use are consistent with the cleared indications.</i>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
		c. A 510(k) number is provided for each component or accessory that received a prior 510(k) clearance. <i>Select "N/A" if the submission states that the component(s)/accessory(ies) does not have a prior 510(k) clearance or the component(s)/accessory(ies) is 510(k) exempt.</i>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
		Comments:			
<b>C.</b>	<b>Substantial Equivalence Discussion</b>				
	14.	Submitter has identified a predicate(s) device	<input checked="" type="checkbox"/>		<input type="checkbox"/>
		a. Predicate's 510(k) number, trade name, and model number (if applicable) provided. For predicates that are preamendments devices, information is provided to document preamendments status. <i>Information regarding <a href="http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidan">documenting preamendment status</a> is available online</i> <i>(<a href="http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidan">http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidan</a>)</i>	<input checked="" type="checkbox"/>		<input type="checkbox"/>

Acceptance Checklist for Traditional 510(k)



*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>					
Submission should be designated RTA if not addressed					
<b>Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.</b>					
			Yes	N/A	No
		<a href="http://www.fda.gov/oc/ComplianceActivities/ucm072746.htm">ce/ComplianceActivities/ucm072746.htm</a>			
	b.	The identified predicate(s) is consistent throughout the submission (i.e., the predicate(s) identified in the Substantial Equivalence section is the same as that listed in the 510(k) Summary (if applicable) and that used in comparative performance testing.	<input checked="" type="checkbox"/>		<input type="checkbox"/>
	Comments: VOL 001 - PAGE 1-12				
	15.	Submission includes a comparison of the following for the predicate(s) and subject device			
	a.	Indications for use	<input checked="" type="checkbox"/>		<input type="checkbox"/>
	b.	Technology, including features, materials, and principles of operation	<input checked="" type="checkbox"/>		<input type="checkbox"/>
	Comments: VOL 001 - PAGE 1-12				
	16.	Submission includes an analysis of why any differences between the subject device and predicate(s) do not render the device NSE (e.g., does not constitute a new intended use; and any differences in technological characteristics are accompanied by information that demonstrates the device is as safe and effective as the predicate and do not raise different questions of safety and effectiveness than the predicate), affect safety or effectiveness, or raise different questions of safety and effectiveness (see section 513(i)(1)(A) of the FD&C Act and 21 CFR 807.87(f)) <i>If there is no difference between the subject and predicate(s) with respect to indications for use or technology, this should be explicitly stated, in which case "N/A" should be selected. Select "No" only if the submission does not include an analysis of differences as described above or a statement that there are no differences. Note that the adequacy of the analysis should be assessed during the substantive review; only the presence of such an analysis is required for acceptance. In addition, note that due to potential differences in</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>				
Submission should be designated RTA if not addressed				
<b>Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.</b>				
	<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>	<b>Yes</b>	<b>N/A</b>	<b>No</b>
	<i>manufacturing that may not be known to the submitter, the fact that no differences are identified does not necessarily mean that no performance testing is needed.</i>			
	Comments:			
<b>D.</b>	<b>Proposed Labeling (see also 21 CFR part 801)</b> <i>If in vitro diagnostic (IVD) device, criteria 17, 18, and 19 may be omitted. These criteria will be omitted from the checklist if "N/A" is selected. IVD labeling is addressed in section 21 below.</i>		<input type="checkbox"/>	
	17. Submission includes proposed package labels and labeling (e.g., instructions for use, package insert, operator's manual) that include a description of the device, its intended use, and the directions for use	<input checked="" type="checkbox"/>		<input type="checkbox"/>
	a. Indications for use are stated in labeling and are identical to Indications for Use form and 510(k) Summary (if 510(k) Summary provided)	<input checked="" type="checkbox"/>		<input type="checkbox"/>
	b. Submission includes directions for use that <ul style="list-style-type: none"> <li>include statements of all conditions, purposes or uses for which the device is intended (e.g., hazards, warnings, precautions, contraindications) (21 CFR 801.5) AND</li> <li>Includes directions for layperson (see 21 CFR 801.5) OR submission states that device qualifies for exemption per 21 CFR 801 Subpart D</li> </ul>	<input checked="" type="checkbox"/>		<input type="checkbox"/>
	Comments: VOL-013			
	18. If indicated for prescription use, labeling includes the prescription use statement (see 21 CFR 801.109(b)(1)) or "Rx only" symbol [See also <a href="#">Alternative to Certain Prescription Device Labeling Requirements</a> ] Select "N/A" if not indicated for prescription use.	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>						
Submission should be designated RTA if not addressed						
<b>Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.</b>						
			Yes	N/A	No	
		<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>				
		Comments: <i>FRONT PAGE IFU - VOL - 013</i>				
	19.	General labeling provisions				
	a.	Labeling includes name and place of business of the manufacturer, packer, or distributor (21 CFR 801.1)	<input checked="" type="checkbox"/>		<input type="checkbox"/>	
	b.	Labeling includes device common or usual name (21 CFR 801.61) <i>Select "N/A" if device is for prescription use only.</i>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
		Comments:				
	20.	a.	If there are requirements regarding labeling, such as special controls, in a device-specific regulation that are applicable to the device, the submission includes labeling to establish that the submitter has followed the device-specific requirement. <i>Select "N/A" if there are no applicable requirements in a device-specific regulation. Select "No" if the submission does not include a rationale for any omitted information. Note that the adequacy of how such requirements have been addressed should be assessed during the substantive review.</i>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
		b.	If there is a device-specific guidance, other than a special controls guidance document, applicable to the device, the submission includes labeling to establish that the submitter has addressed the recommendations or otherwise has met the applicable statutory or regulatory criteria through an alternative approach. <i>Select "N/A" if there is no applicable device-specific guidance. Select "No" if the submission does not include a rationale for any omitted information or any alternative approach as outlined above. Note that the adequacy of how recommendations in a device-specific guidance have been addressed should be assessed during the substantive review.</i>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>					
Submission should be designated RTA if not addressed					
<b>Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.</b>					
			Yes	N/A	No
<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>					
	c.	<p>If there is a special controls document applicable to the device, the submission includes labeling to establish that the submitter has complied with the particular mitigation measures set forth in the special controls document or uses alternative mitigation measures but provides a rationale to demonstrate that those alternative measures identified by the firm will provide at least an equivalent assurance of safety and effectiveness.</p> <p><i>Select "N/A" if there is no applicable special controls document. Select "No" if the submission does not include a rationale for any omitted information or any alternative approach as outlined above. Note that the adequacy of how mitigation measures in a special controls document have been addressed should be assessed during the substantive review.</i></p>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
		Comments:			
	21.	<p>If the device is an in vitro diagnostic device, provided labeling includes all applicable information required per 21 CFR 809.10.</p> <p><i>Select "N/A" if not an in vitro diagnostic device.</i></p>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
<b>E.</b>	<p><b>Sterilization</b></p> <p><i>If in vitro diagnostic (IVD) device and sterilization is not applicable, select "N/A." The criteria in this section will be omitted from the checklist if "N/A" is selected.</i></p>			<input checked="" type="checkbox"/>	
	<p>Submission states that the device and/or accessories are: <i>(one of the below must be checked)</i></p> <p><input type="checkbox"/> provided sterile</p> <p><input type="checkbox"/> provided non-sterile but sterilized by the end user</p> <p><input type="checkbox"/> non-sterile when used</p> <p>This information will determine whether and what type of additional information may be necessary for a substantial equivalence determination.</p> <p><i>If "non-sterile when used" is selected, the sterility-related criteria below are omitted from</i></p>				<input type="checkbox"/>

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>					
Submission should be designated RTA if not addressed					
<b>Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.</b>					
	<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>	<b>Yes</b>	<b>N/A</b>	<b>No</b>	
	<i>the checklist.</i> <i>If information regarding the sterility status of the device is not provided, select "No."</i>				
	Comments: <b>STANDALONE SOFTWARE</b>				
	22.	Assessment of the need for sterilization information			
	a.	Identification of device, and/or accessories, and/or components that are provided sterile.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	b.	Identification of device, and/or accessories, and/or components that are end user sterilized	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	c.	Identification of device, and/or accessories, and/or components that are reusable and cleaning/disinfection instructions are provided.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Comments:			
	23.	If the device, and/or accessory, and/or a component is provided sterile: <i>Select "N/A" if no part of the device, accessories, or components is provided sterile, otherwise complete a-e below.</i>		<input type="checkbox"/>	
	a.	Sterilization method is stated for each component (including parameters such as dry time for steam sterilization, radiation dose, etc.)	<input type="checkbox"/>		<input type="checkbox"/>
	b.	A description of method to validate the sterilization parameters (e.g., half-cycle method and full citation of FDA-recognized standard, including date) is provided for each proposed sterilization method. <i>Note, the sterilization validation report is not required.</i>	<input type="checkbox"/>		<input type="checkbox"/>
	c.	For devices sterilized using chemical sterilants such as ethylene oxide (EO) and hydrogen peroxide, submission states maximum	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>						
Submission should be designated RTA if not addressed						
<b>Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.</b>						
				Yes	N/A	No
		<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>				
		levels of sterilant residuals remaining on the device and sterilant residual limits. <i>Select "N/A" if not sterilized using chemical sterilants.</i>				
	d.	Submission includes description of packaging and packaging contents (e.g., if multiple devices are included within the same package, Tyvek packaging, etc.)		<input type="checkbox"/>		<input type="checkbox"/>
	e.	Sterility Assurance Level (SAL) stated		<input type="checkbox"/>		<input type="checkbox"/>
		Comments:				
	24.	If the device, and/or accessory, and/or a component is end user sterilized: <i>Select "N/A" if no part of the device, accessories, or components are end user sterilized, otherwise complete a-d below.</i>			<input type="checkbox"/>	
	a.	Sterilization method is stated for each component (including parameters such as dry time for steam sterilization, radiation dose, etc.)		<input type="checkbox"/>		<input type="checkbox"/>
	b.	A description of method to validate the sterilization parameters (e.g., half-cycle method and full citation of FDA-recognized standard, including date) is provided for each proposed sterilization method. <i>Note, the sterilization validation is not required.</i>		<input type="checkbox"/>		<input type="checkbox"/>
	c.	Submission includes description of packaging and packaging contents (e.g., if multiple devices are included within the same package, Tyvek packaging, etc.)		<input type="checkbox"/>		<input type="checkbox"/>
	d.	Submission includes sterilization instructions for end user		<input type="checkbox"/>		<input type="checkbox"/>
		Comments:				
	25.	a. If there are requirements regarding sterility, such as special		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>							
Submission should be designated RTA if not addressed							
<b>Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.</b>							
<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>					Yes	N/A	No
			controls, in a device-specific regulation that are applicable to the device, the submission includes sterility information to establish that the submitter has followed the device-specific requirement. <i>Select "N/A" if there are no applicable requirements in a device-specific regulation. Select "No" if the submission does not include a rationale for any omitted information. Note that the adequacy of how such requirements have been addressed should be assessed during the substantive review.</i>				
		b.	If there is a device-specific guidance, other than a special controls guidance document, applicable to the device, the submission includes sterility information to establish that the submitter has addressed the recommendations or otherwise has met the applicable statutory or regulatory criteria through an alternative approach. <i>Select "N/A" if there is no applicable device-specific guidance. Select "No" if the submission does not include a rationale for any omitted information or any alternative approach as outlined above. Note that the adequacy of how recommendations in a device-specific guidance have been addressed should be assessed during the substantive review.</i>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		c.	If there is a special controls document applicable to the device, the submission includes sterility information to establish that the submitter has complied with the particular mitigation measures set forth in the special controls document or uses alternative mitigation measures but provides a rationale to demonstrate that those alternative measures identified by the firm will provide at least an equivalent assurance of safety and effectiveness. <i>Select "N/A" if there is no applicable special controls document. Select "No" if the submission does not include a rationale for any omitted information or any alternative approach as outlined above. Note that the adequacy of how mitigation measures in a</i>				

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>					
Submission should be designated RTA if not addressed					
<b>Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.</b>					
			Yes	N/A	No
		<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>			
		<i>special controls document have been addressed should be assessed during the substantive review.</i>			
		Comments:			
<b>F.</b>	<b>Shelf Life</b>				
	26.	Proposed shelf life/ expiration date stated <i>Select "N/A" if the device is not provided sterile and the submitter states that storage conditions could not affect device safety or effectiveness.</i>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
		Comments:			
	27.	For sterile device, submission includes summary of methods used to establish that device sterility will remain substantially equivalent to that of the predicate through the proposed shelf life, or a rationale for why testing to establish shelf life is not applicable. <i>Select "N/A" if the device is not provided sterile.</i>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
		Comments:			
	28.	Submission includes summary of methods used to establish that device performance is not adversely affected by aging and therefore device performance will remain substantially equivalent to that of the predicate, or includes a rationale for why the storage conditions are not expected to affect device safety or effectiveness.	<input type="checkbox"/>		<input type="checkbox"/>
		Comments:			
<b>G.</b>	<b>Biocompatibility</b> <i>If in vitro diagnostic (IVD) device, select "N/A." The criteria in this section will be omitted from the checklist if "N/A" is selected.</i>				<input checked="" type="checkbox"/>
	Submission states that there: <i>(one of the below must be checked)</i>				<input type="checkbox"/>
	<input type="checkbox"/> are				

Acceptance Checklist for Traditional 510(k)



*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>				
Submission should be designated RTA if not addressed				
<b>Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.</b>				
		Yes	N/A	No
	<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>			
	<input checked="" type="checkbox"/> are not direct or indirect (e.g., through fluid infusion) patient-contacting components.  This information will determine whether and what type of additional information may be necessary for a substantial equivalence determination. <i>If "are not" is selected, the biocompatibility-related criteria below are omitted from the checklist. If information regarding whether the device is patient-contacting is not provided, select "No."</i>			
	Comments: VOL_001 PAGE 1-18			
	29. Submission includes list of patient-contacting device components and associated materials of construction, including identification of color additives, if present	<input type="checkbox"/>		<input type="checkbox"/>
	Comments:			
	30. Submission identifies contact classification (e.g., surface-contacting, less than 24 hour duration)	<input type="checkbox"/>		<input type="checkbox"/>
	Comments:			
	31. Biocompatibility assessment of patient-contacting components  Submission includes: Test protocol (including identification and description of test article), methods, pass/fail criteria, and results provided for each completed test, OR a statement that biocompatibility testing is not needed with a rationale (e.g., materials and manufacturing/processing are identical to the predicate).	<input type="checkbox"/>		<input type="checkbox"/>
	Comments:			
<b>H.</b>	<b>Software</b>			

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>				
Submission should be designated RTA if not addressed				
<b>Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.</b>				
		Yes	N/A	No
	<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>			
	Submission states that the device: <i>(one of the below must be checked)</i> <input checked="" type="checkbox"/> does <input type="checkbox"/> does not contain software/firmware.  This information will determine whether and what type of additional information may be necessary for a substantial equivalence determination. <i>If "does not" is selected, the software-related criterion is omitted from the checklist. If information regarding whether the device contains software is not provided, select "No."</i>			<input type="checkbox"/>
	Comments:			
	32. Submission includes a statement of software level of concern and rationale for the software level of concern	<input checked="" type="checkbox"/>		<input type="checkbox"/>
	Comments: VOL 002			
	33. All applicable software documentation provided based on level of concern identified by the submitter, as described in <a href="#">Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices</a> , or the submission includes information to establish that the submitter has otherwise met the applicable statutory or regulatory criteria through an alternative approach (i.e., the submitter has identified an alternate approach with a rationale).	<input checked="" type="checkbox"/>		<input type="checkbox"/>
	Comments: VOL 003 To Vol 013			
<b>I.</b>	<b>EMC and Electrical Safety</b>			
	Submission states that the device: <i>(one of the below must be checked)</i> <input type="checkbox"/> does <input checked="" type="checkbox"/> does not require EMC and Electrical Safety evaluation.			<input type="checkbox"/>

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>				
Submission should be designated RTA if not addressed				
<b>Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.</b>				
		Yes	N/A	No
	<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>			
	This information will determine whether and what type of additional information may be necessary for a substantial equivalence determination. <i>If "does not" is selected, the EMC-related and Electrical Safety-related criteria below are omitted from the checklist. If information regarding whether the device requires EMC and Electrical Safety evaluation is not provided, select "No."</i>			
	Comments: VOL - 001 PAGE 1-16			
34.	Submission includes evaluation of electrical safety (e.g., per IEC 60601-1, or equivalent FDA-recognized standard, and if applicable, the device-specific standard), OR submission includes electrical safety evaluation using methods or standards that are not FDA-recognized and submission includes information to establish that the submitter has otherwise met the applicable statutory or regulatory criteria through this alternative approach (i.e., the submitter has identified alternate methods or standards with a rationale).	<input type="checkbox"/>		<input type="checkbox"/>
	Comments:			
35.	Submission includes evaluation of electromagnetic compatibility (e.g., per IEC 60601-1-2 or equivalent FDA-recognized standard and if applicable, the device-specific standard) OR submission includes electromagnetic compatibility evaluation using methods or standards that are not FDA-recognized and submission includes information to establish that the submitter has otherwise met the applicable statutory or regulatory criteria through this alternative approach (i.e., the submitter has identified alternate methods or standards with a rationale).	<input type="checkbox"/>		<input type="checkbox"/>
	Comments:			

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>					
Submission should be designated RTA if not addressed					
<b>Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.</b>					
			Yes	N/A	No
	<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>				
<b>J.</b>	<b>Performance Data – General</b> <i>If in vitro diagnostic (IVD) device, select "N/A." The criteria in this section will be omitted from the checklist if "N/A" is selected. Performance data criteria relating to IVD devices will be addressed in Section K.</i>			<input type="checkbox"/>	
	Comments:				
	36. Full test report is provided for each completed test. A full test report includes: objective of the test, description of the test methods and procedures, study endpoint(s), pre- defined pass/fail criteria, results summary, conclusions, and an explanation of how the data generated from the test supports a finding of substantial equivalence. <i>Full test reports provided for all completed tests/evaluations (e.g., bench evaluations, comparative performance tests, etc.). Select "N/A" if the submission does not include performance data.</i>		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Comments: <i>VAL OLO</i>				
	37. a. If there are requirements regarding performance data, such as special controls, in a device-specific regulation that are applicable to the device, the submission includes performance data to establish that the submitter has followed the device-specific requirement. <i>Select "N/A" if there are no applicable requirements in a device-specific regulation. Select "No" if the submission does not include a rationale for any omitted information. Note that the adequacy of how such requirements have been addressed should be assessed during the substantive review.</i>		<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
	b. If there is a device-specific guidance, other than a special controls guidance document, applicable to the device, the submission includes performance data to establish that the submitter has addressed the recommendations or otherwise has met the applicable statutory or regulatory criteria through an alternative		<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>						
Submission should be designated RTA if not addressed						
<b>Check “Yes” if item is present, “N/A” if it is not needed and “No” if it is not included but needed.</b>						
				Yes	N/A	No
<ul style="list-style-type: none"> <li>Any “No” answer will result in a “Refuse to Accept” decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>						
			approach. <i>Select “N/A” if there is no applicable device-specific guidance. Select “No” if the submission does not include a rationale for any omitted information or any alternative approach as outlined above. Note that the adequacy of how recommendations in a device-specific guidance have been addressed should be assessed during the substantive review.</i>			
		c.	If there is a special controls document applicable to the device, the submission includes performance data to establish that the submitter has complied with the particular mitigation measures set forth in the special controls document or uses alternative mitigation measures but provides a rationale to demonstrate that those alternative measures identified by the firm will provide at least an equivalent assurance of safety and effectiveness. <i>Select “N/A” if there is no applicable special controls document. Select “No” if the submission does not include a rationale for any omitted information or any alternative approach as outlined above. Note that the adequacy of how mitigation measures in a special controls document have been addressed should be assessed during the substantive review.</i>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Comments:						
	38.		If literature is referenced in the submission, submission includes: <i>Select “N/A” if the submission does not reference literature. Note that the applicability of the referenced article to support a substantial equivalence finding should be assessed during the substantive review; only the presence of a discussion is required to support acceptance.</i>		<input checked="" type="checkbox"/>	
		a.	Legible reprints or a summary of each article	<input type="checkbox"/>		<input type="checkbox"/>
		b.	Discussion of how each article is applicable to support the	<input type="checkbox"/>		<input type="checkbox"/>

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>				
Submission should be designated RTA if not addressed				
<b>Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.</b>				
	<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>	<b>Yes</b>	<b>N/A</b>	<b>No</b>
	substantial equivalence of the subject device to the predicate.			
	Comments:			
	39. For each completed nonclinical (i.e., animal) study conducted, <i>Select "N/A" if no animal study was conducted. Note that this section does not address biocompatibility evaluations, which are assessed in Section G of the checklist,</i>		<input checked="" type="checkbox"/>	
	a. Submission includes a study protocol which includes all elements as outlined in 21 CFR 58.120	<input type="checkbox"/>		<input type="checkbox"/>
	b. Submission includes final study report which includes all elements outlined in 21 CFR 58.185	<input type="checkbox"/>		<input type="checkbox"/>
	c. Submission contains a statement that the study was conducted in compliance with applicable requirements in the GLP regulation (21 CFR Part 58), or, if the study was not conducted in compliance with the GLP regulation, the submission explains why the noncompliance would not impact the validity of the study data provided to support a substantial equivalence determination.	<input type="checkbox"/>		<input type="checkbox"/>
	Comments:			
<b>K.</b>	<b>Performance Characteristics – In Vitro Diagnostic Devices Only (see also 21 CFR 809.10(b)(12))</b>			
	Submission indicates that device: <i>(one of the below must be checked)</i> <input type="checkbox"/> is <input checked="" type="checkbox"/> is not an in vitro diagnostic device (IVD). <i>If "is not" is selected, the performance data-related criteria below are omitted from the checklist.</i>			
	Comments:			

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>					
Submission should be designated RTA if not addressed					
Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.					
		<ul style="list-style-type: none"> <li>Any "No" answer will result in a "Refuse to Accept" decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>	Yes	N/A	No
	40.	Submission includes the following studies, as appropriate for the device type, including associated protocol descriptions, study results and line data:			
	a.	Precision/reproducibility	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	b.	Accuracy (includes as appropriate linearity; calibrator or assay traceability; calibrator and/or assay stability protocol and acceptance criteria; assay cut-off; method comparison or comparison to clinical outcome; matrix comparison; and clinical reference range or cutoff).	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	c.	Sensitivity (detection limits, LoB, LoD, LoQ where relevant for the device type).	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	d.	Analytical specificity	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Comments:			
	41.	a. If there are requirements regarding performance data, such as special controls, in a device-specific regulation that are applicable to the device, the submission includes performance data to establish that the submitter has followed the device-specific requirement. <i>Select "N/A" if there are no applicable requirements in a device-specific regulation. Select "No" if the submission does not include a rationale for any omitted information. Note that the adequacy of how such requirements have been addressed should be assessed during the substantive review.</i>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		b. If there is a device-specific guidance, other than a special controls guidance document, applicable to the device, the submission includes performance data to establish that the submitter has addressed the recommendations or otherwise has met the	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Acceptance Checklist for Traditional 510(k)

*Contains Nonbinding Recommendations*

<b><u>Elements of a Complete Submission (RTA Items)</u></b> <b><u>(21 CFR 807.87 unless otherwise indicated)</u></b>							
Submission should be designated RTA if not addressed							
<b>Check “Yes” if item is present, “N/A” if it is not needed and “No” if it is not included but needed.</b>							
<ul style="list-style-type: none"> <li>Any “No” answer will result in a “Refuse to Accept” decision.</li> <li>Each element on the checklist should be addressed within the submission. The submitter may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criterion is considered present (Yes). An assessment of the rationale will be considered during the review of the submission.</li> </ul>					Yes	N/A	No
			applicable statutory or regulatory criteria through an alternative approach. <i>Select “N/A” if there is no applicable device-specific guidance. Select “No” if the submission does not include a rationale for any omitted information or any alternative approach as outlined above. Note that the adequacy of how recommendations in a device-specific guidance have been addressed should be assessed during the substantive review.</i>				
		c.	If there is a special controls document applicable to the device, the submission includes performance data to establish that the submitter has complied with the particular mitigation measures set forth in the special controls document or uses alternative mitigation measures but provides a rationale to demonstrate that those alternative measures identified by the firm will provide at least an equivalent assurance of safety and effectiveness. <i>Select “N/A” if there is no applicable special controls document. Select “No” if the submission does not include a rationale for any omitted information or any alternative approach as outlined above. Note that the adequacy of how mitigation measures in a special controls document have been addressed should be assessed during the substantive review.</i>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Comments:							

Acceptance Checklist for Traditional 510(k)



Department of Health and Human Services  
Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
(To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional                       Special                       Abbreviated

STANDARD TITLE <sup>1</sup>

IEC 62304 Medical Device Software - Software Life Cycle Processes, 08/20/2012

**Please answer the following questions**

Yes      No

Is this standard recognized by FDA <sup>2</sup>? .....      

FDA Recognition number <sup>3</sup> ..... # 13-8

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report <sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) <sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
If yes, report these exclusions in the summary report table.

Is there an FDA guidance <sup>6</sup> that is associated with this standard? .....         
If yes, was the guidance document followed in preparation of this 510k? .....      

Title of guidance: General Principles of Software Validation; Final Guidance for Industry and FDA Staff

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/Standards/default.htm>

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and

address of the test laboratory or certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/default.htm>

EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE		
STANDARD TITLE IEC 62304 Medical Device Software - Software Life Cycle Processes, 08/20/2012		
CONFORMANCE WITH STANDARD SECTIONS*		
SECTION NUMBER 4	SECTION TITLE Quality Management System	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
SECTION NUMBER 5	SECTION TITLE Software development process	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
SECTION NUMBER 6	SECTION TITLE Software maintenance process	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
<p>* For completeness list all sections of the standard and indicate whether conformance is met. If a section is not applicable (N/A) an explanation is needed under "justification." Some standards include options, so similar to deviations, the option chosen needs to be described and adequately justified as appropriate for the subject device. Explanation of all deviations or description of options selected when following a standard is required under "type of deviation or option selected," "description" and "justification" on the report. More than one page may be necessary.</p> <p>♦ Types of deviations can include an exclusion of a section in the standard, a deviation brought out by the FDA supplemental information sheet (SIS), a deviation to adapt the standard to the device, or any adaptation of a section.</p>		
<p>This section applies only to requirements of the Paperwork Reduction Act of 1995.</p> <p><b>*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.*</b></p> <p>The burden time for this collection of information is estimated to average 1 hour per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:</p> <div style="display: flex; justify-content: space-between;"> <div style="width: 45%;"> <p>Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff <a href="mailto:PRAStaff@fda.hhs.gov">PRAStaff@fda.hhs.gov</a></p> </div> <div style="width: 45%; text-align: right;"> <p><i>"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number."</i></p> </div> </div>		

EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE		
STANDARD TITLE IEC 62304 Medical Device Software - Software Life Cycle Processes, 08/20/2012		
CONFORMANCE WITH STANDARD SECTIONS*		
SECTION NUMBER 7	SECTION TITLE Software risk management process	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
SECTION NUMBER 8	SECTION TITLE Software configuration management process	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
SECTION NUMBER 9	SECTION TITLE Software problem resolution process	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
SECTION NUMBER	SECTION TITLE	CONFORMANCE?
		<input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
<p>* For completeness list all sections of the standard and indicate whether conformance is met. If a section is not applicable (N/A) an explanation is needed under "justification." Some standards include options, so similar to deviations, the option chosen needs to be described and adequately justified as appropriate for the subject device. Explanation of all deviations or description of options selected when following a standard is required under "type of deviation or option selected," "description" and "justification" on the report. More than one page may be necessary.</p> <p>♦ Types of deviations can include an exclusion of a section in the standard, a deviation brought out by the FDA supplemental information sheet (SIS), a deviation to adapt the standard to the device, or any adaptation of a section.</p>		

Department of Health and Human Services  
Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
(To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional                       Special                       Abbreviated

STANDARD TITLE <sup>1</sup>

ISO 14971 Medical devices - Application of Risk Management To Medical Devices 16/05/2012

**Please answer the following questions**

Yes      No

Is this standard recognized by FDA <sup>2</sup>? .....      

FDA Recognition number <sup>3</sup> ..... #5-40

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report <sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....         
If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) <sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....         
If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
If yes, report these exclusions in the summary report table.

Is there an FDA guidance <sup>6</sup> that is associated with this standard? .....         
If yes, was the guidance document followed in preparation of this 510k? .....      

Title of guidance: \_\_\_\_\_

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/Standards/default.htm>

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and

address of the test laboratory or certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/default.htm>

<b>EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE</b>		
STANDARD TITLE ISO 14971 Medical devices - Application of Risk Management To Medical Devices 16/05/2012		
<b>CONFORMANCE WITH STANDARD SECTIONS*</b>		
SECTION NUMBER 3	SECTION TITLE General requirements for risk management	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
SECTION NUMBER 4	SECTION TITLE Risk analysis	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
SECTION NUMBER 5	SECTION TITLE Risk evaluation	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
<p>* For completeness list all sections of the standard and indicate whether conformance is met. If a section is not applicable (N/A) an explanation is needed under "justification." Some standards include options, so similar to deviations, the option chosen needs to be described and adequately justified as appropriate for the subject device. Explanation of all deviations or description of options selected when following a standard is required under "type of deviation or option selected," "description" and "justification" on the report. More than one page may be necessary.</p> <p>♦ Types of deviations can include an exclusion of a section in the standard, a deviation brought out by the FDA supplemental information sheet (SIS), a deviation to adapt the standard to the device, or any adaptation of a section.</p>		
<p>This section applies only to requirements of the Paperwork Reduction Act of 1995.</p> <p><b>*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.*</b></p> <p>The burden time for this collection of information is estimated to average 1 hour per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:</p> <p style="margin-left: 40px;">Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff <a href="mailto:PRAStaff@fda.hhs.gov">PRAStaff@fda.hhs.gov</a></p> <p style="margin-left: 40px;"><i>"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number."</i></p>		

EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE		
STANDARD TITLE ISO 14971 Medical devices - Application of Risk Management To Medical Devices 16/05/2012		
CONFORMANCE WITH STANDARD SECTIONS*		
SECTION NUMBER 6	SECTION TITLE Risk control	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
SECTION NUMBER 7	SECTION TITLE Evaluation of overall residual risk acceptability	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
SECTION NUMBER 8	SECTION TITLE Risk management report	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
SECTION NUMBER 9	SECTION TITLE Production and post-production information	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
<p>* For completeness list all sections of the standard and indicate whether conformance is met. If a section is not applicable (N/A) an explanation is needed under "justification." Some standards include options, so similar to deviations, the option chosen needs to be described and adequately justified as appropriate for the subject device. Explanation of all deviations or description of options selected when following a standard is required under "type of deviation or option selected," "description" and "justification" on the report. More than one page may be necessary.</p> <p>♦ Types of deviations can include an exclusion of a section in the standard, a deviation brought out by the FDA supplemental information sheet (SIS), a deviation to adapt the standard to the device, or any adaptation of a section.</p>		

Department of Health and Human Services  
Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
(To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional                       Special                       Abbreviated

STANDARD TITLE <sup>1</sup>

ISO 13485 Medical devices - Quality management systems - Requirements for regulatory proces, 2003-07-15

**Please answer the following questions**

Yes      No

Is this standard recognized by FDA <sup>2</sup>? .....      

FDA Recognition number <sup>3</sup> ..... # \_\_\_\_\_

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report <sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard?.....         
If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) <sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS?.....         
If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
If yes, report these exclusions in the summary report table.

Is there an FDA guidance <sup>6</sup> that is associated with this standard?.....         
If yes, was the guidance document followed in preparation of this 510k? .....      

Title of guidance: \_\_\_\_\_

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/Standards/default.htm>

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and

address of the test laboratory or certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/default.htm>

**EXTENT OF STANDARD CONFORMANCE  
SUMMARY REPORT TABLE**

STANDARD TITLE  
ISO 13485 Medical devices - Quality management systems - Requirements for regulatory proces, 2003-07-15

**CONFORMANCE WITH STANDARD SECTIONS\***

SECTION NUMBER 4	SECTION TITLE Quality management system	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
---------------------	--	--

TYPE OF DEVIATION OR OPTION SELECTED ♦

DESCRIPTION

JUSTIFICATION

SECTION NUMBER 5	SECTION TITLE Management responsibility	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
---------------------	--	--

TYPE OF DEVIATION OR OPTION SELECTED ♦

DESCRIPTION

JUSTIFICATION

SECTION NUMBER 6	SECTION TITLE Resource management	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
---------------------	--------------------------------------	--

TYPE OF DEVIATION OR OPTION SELECTED ♦

DESCRIPTION

JUSTIFICATION

\* For completeness list all sections of the standard and indicate whether conformance is met. If a section is not applicable (N/A) an explanation is needed under "justification." Some standards include options, so similar to deviations, the option chosen needs to be described and adequately justified as appropriate for the subject device. Explanation of all deviations or description of options selected when following a standard is required under "type of deviation or option selected," "description" and "justification" on the report. More than one page may be necessary.

♦ Types of deviations can include an exclusion of a section in the standard, a deviation brought out by the FDA supplemental information sheet (SIS), a deviation to adapt the standard to the device, or any adaptation of a section.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 1 hour per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRAStaff@fda.hhs.gov](mailto:PRAStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number."*



EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE		
STANDARD TITLE ISO 13485 Medical devices - Quality management systems - Requirements for regulatory proces, 2003-07-15		
<b>CONFORMANCE WITH STANDARD SECTIONS*</b>		
SECTION NUMBER 7	SECTION TITLE Product realization	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
SECTION NUMBER 8	SECTION TITLE Measurement, analysis and improvement	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
SECTION NUMBER	SECTION TITLE	CONFORMANCE? <input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
SECTION NUMBER	SECTION TITLE	CONFORMANCE? <input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
TYPE OF DEVIATION OR OPTION SELECTED ♦		
DESCRIPTION		
JUSTIFICATION		
<p>* For completeness list all sections of the standard and indicate whether conformance is met. If a section is not applicable (N/A) an explanation is needed under "justification." Some standards include options, so similar to deviations, the option chosen needs to be described and adequately justified as appropriate for the subject device. Explanation of all deviations or description of options selected when following a standard is required under "type of deviation or option selected," "description" and "justification" on the report. More than one page may be necessary.</p> <p>♦ Types of deviations can include an exclusion of a section in the standard, a deviation brought out by the FDA supplemental information sheet (SIS), a deviation to adapt the standard to the device, or any adaptation of a section.</p>		

## Level of concern

Level of Concern is determined by answering the key questions listed in [Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices](#), Table 1 and Table 2.

**Table 1 Major Level of Concern**

If the answer to any one question below is Yes, the Level of Concern for the Software Device is likely to be Major.
Does the Software Device qualify as Blood Establishment Computer Software? <b>No</b>
(Blood Establishment Computer Software is defined as software products intended for use in the manufacture of blood and blood components or for the maintenance of data that blood establishment personnel use in making decisions regarding the suitability of donors and the release of blood or blood components for transfusion or further manufacture.)
Is the Software Device intended to be used in combination with a drug or biologic? <b>No</b>
Is the Software Device an accessory to a medical device that has a Major Level of Concern? <b>No</b>
Prior to mitigation of hazards, could a failure of the Software Device result in death or serious injury, either to a patient or to a user of the device? Examples of this include the following:
Does the Software Device control a life supporting or life sustaining function? <b>No</b>
Does the Software Device control the delivery of potentially harmful energy that could result in death or serious injury, such as radiation treatment systems, defibrillators, and ablation generators? <b>No</b>
Does the Software Device control the delivery of treatment or therapy such that an error or malfunction could result in death or serious injury? <b>No</b>
Does the Software Device provide diagnostic information that directly drives a decision regarding treatment or therapy, such that if misapplied it could result in serious injury or death? <b>No</b>
Does the Software Device provide vital signs monitoring and alarms for potentially life threatening situations in which medical intervention is necessary? <b>No</b>

Based on the above Dental application software is not Major Level of Concern.



3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

**Table 2 Moderate Level of Concern**

If the Software Device is not Major Level of Concern and the answer to any one question below is Yes, the Level of Concern is likely to be Moderate.
Is the Software Device an accessory to a medical device that has a Moderate Level of Concern? <b>No</b>
Prior to mitigation of hazards, could a failure of the Software Device result in Minor Injury, either to a patient or to a user of the device? <b>Yes</b>
Could a malfunction of, or a latent design flaw in, the Software Device lead to an erroneous diagnosis or a delay in delivery of appropriate medical care that would likely lead to Minor Injury? <b>No</b>

The answers to all of the questions in Table 1 are **No**, which is why Level of Concern for the 3shape Abutment Designer™ Software is determined to be **Moderate**.

---

## Software Description

### 1. Introduction

The 3Shape Abutment Designer™ Software is intended for designing endosseous dental implant abutments based on digital scans of 3D models.

The output of the device is a computer file containing the abutment(s) in digital form which can be used by 3<sup>rd</sup> party providers of milling machines and/or rapid prototyping machines to manufacture physical abutments. The 3Shape Abutment Designer™ Software does not provide any means of creating physical abutments, nor does it consist of any physical parts.

The 3Shape Abutment Designer™ Software can be run on properly configured "off the shelf" PC hardware running Microsoft® Windows and using a standard keyboard and mouse.

The Software is implemented in the Delphi programming language.

For summary of functional requirements of the Software, see "VOL 005 Software Requirement Specification" of this submission.

**Note:** 3Shape Abutment Designer™ is an add on module to 3Shape Dental System™ which is classified as 872.3661 Optical Impression Systems for CAD/CAM, Product Code NOF (510(K) Exempt). Throughout this submission there will be references to 3Shape Dental System™, but only the documentation relevant to 3Shape Abutment Designer™ has been included.



3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

## **2. Intended Use and Operational Environment**

The 3Shape Abutment Designer™ Software is intended for designing endosseous dental implant abutments based on digital scans of 3D models.

Intended users are Dental Practitioners and Dental Laboratory staff.

Intended Operational Environment is Dental Laboratories.

## 3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

**3. Definitions**


There are four types of Abutments supported by the 3Shape Dental System: Customized, Anatomical, Screw Retained Crown, and Wax up.

3Shape uses the term "Abutment" to denote all four types. Digitally created abutments are classified as an Endosseous dental implant abutment.

The 3Shape Abutment Designer™ only supports modification of the Abutment itself as exemplified on *Figure 1: Customized Abutment parts*.

A digital representation of the Implant is required for the software to work, but only to ensure the Abutment will fit on the Seating Geometry / Interface of the Implant. The 3Shape Abutment Designer™ does not provide any means to design, alter, or manufacture any part of the Implant including, but not limited to, the abutment to implant interface.

**Figure 1: Customized Abutment parts**

**Terms used**

The 3Shape Abutment Designer™ Software supports four types of Abutments however this only has impact on the initial guess of abutment shape as well as slight differences in the User Interface for designing the Abutment.



3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

There is no difference between the four types on abutment function, manufacturing requirements, output types, etc.

The output of any abutment designed is a computer file of the surface model.

Please find a detailed description in *Table 2: Abutment Types*

**Table 2: Abutment Types**

Type	Description
Customized Abutment	<p>The Customized Abutment comes in 3 subtypes: Custom, Robotic, and Bar Interface.</p> <p>All provide the user with the ability to free hand modify the Top Cap and Emergence Profile as depicted on <i>Figure 1: Customized Abutment parts</i>.</p> <p>The Custom Abutment initial shape guess is based on an anatomic heuristic</p>  <p>The Robotic initial shape guess is based on a classic standard abutment.</p>  <p>The Bar Interface initial shape guess is based on a cylinder.</p>

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

	
<p>Anatomical Abutment</p>	<p>This type of abutment comes with an initial shape guess that matches a crown.</p> 
<p>Screw Retained Crown</p>	<p>This type of abutment comes with an initial shape guess that can be used as a crown and has a screw hole in the centre.</p> 
<p>Wax up abutment</p>	<p>This type of abutment comes with an initial shape guess based on a scan of a wax model. Otherwise, it is equivalent to Anatomical abutment.</p>

Also, throughout the submission the following terms are applicable.

**Table 1: Terms**

Term	Definition
Provider	Any creator of an Implant Library.
Implant Library	A digital representation of an Implant System
Implant System	A physical Dental Implant compatible with one of more abutments
Abutment	Endosseous dental implant abutment as regulated under 872.3630





3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

Milling Center	<p>An industry term for a large scale provider of Milling and/or 3D Printing. Typically, the actual manufacturing is geographically centralized, which makes digital transfer of CAD designed surface files convenient.</p> <p>3Shape does not own or operate any milling centers.</p>
Local Milling	<p>When the Dental Lab and milling machinery is co located.</p> <p>3Shape does not manufacture or market milling machinery.</p>



3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

#### 4. Ensuring Regulatory Compliance of Abutments

Abutments are designed to interface with Dental Implants and 3Shape Abutment Designer™ requires a digital representation of the Dental Implant in order to create a functional Abutment.


The digital representations of Dental Implants (called Implant Systems) is collected in files called Implant Libraries.

In the U.S. only Abutments designed against Dental Implants with a 510(k) clearance are allowed in the 3Shape Abutment Designer™ (see *Design stops in the software* below).

Providers of Implant Libraries to be used in the U.S. must supply 3Shape with written documentation in order to have their Libraries activated in the software.

It is the responsibility of the Provider to ensure that the digital representation correctly represent the cleared physical parts.

#### Design stops in the software





3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the section header and ending above the footer. The text "(b)(4)" is written in red at the top left corner of this redacted area.



3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

---

## Device Hazard Analysis


(b)(4)



**PRODUCT RISK ASSESSMENT**

(b)(4)





(b)(4)

(b)(4)



## Software Requirement Specification

### 1. Introduction

Please refer to the Software Requirement Specifications enclosed in this volume.

Please note that the Validation Protocol templates are included in the Requirement Specifications.

For links between the Requirements Specifications and the Software Specifications, please see "VOL 008 Traceability Analysis" enclosed in this submission.

**Note:** Some of these documents may have requirements referring to other areas than 3Shape Abutment Designer™. See the Trace matrix for reference to which requirements are relevant for 3Shape Abutment Designer™.

### 2. Table of Contents

(b)(4)	Configuration
	Sending and receiving
	Manufacturing Interfaces
	User Manual
	Material import export
	Parametric Abutment Design
	Wax-up Design
	Implant Bars



Requirement Specification - (b)(4) 

(b)(4)



(b)(4)



(b)(4)





Requirement Specification - (b)(4)

(b)(4)



(b)(4)





Requirement Specification - (b)(4)

(b)(4)

A large black rectangular redaction box covers the entire central portion of the page, obscuring all text and graphics that would otherwise be present.

(b)(4)



(b)(4)







Requirement Specification - (b)(4)

(b)(4)



(b)(4)





Requirement Specification - (b)(4)

(b)(4)



(b)(4)





Requirement Specification - (b)(4)

(b)(4)




(b)(4)



(b)(4)







(b)(4)




(b)(4)



Requirement Specification - (b)(4) 


(b)(4)





Requirement Specification - (b)(4)

(b)(4)



(b)(4)



(b)(4)



(b)(4)





(b)(4)




---

## Architecture Design Chart

(b)(4)





(b)(4)



(b)(4)




TITLE: ARCHITECTURE DESIGN CHART

(b)(4)



TITLE: ARCHITECTURE DESIGN CHART



## Software Design Specification

### 1. Introduction

Please refer to the Software Design Specification enclosed in this volume.

Please note that the Verification Protocol templates are included in the Software Specifications.

For links between the Requirements Specifications and the Software Specifications, please see "VOL 008 Traceability Analysis" enclosed in this submission.


### 2. Software Design Specification, Table of Contents

(b)(4)	Sending and receiving
	Manufacturing Interfaces
	User Manual
	Material Import and Export
	Abutment Design Configuration
	Parametric Abutment Design
	Wax up Design
	Implant Bar Design
	RawSTL Manufacturing Output



Software Specification – (b)(4)

(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)





Software Specification – (b)(4)

(b)(4)


A large black rectangular redaction box covers the entire central portion of the page, obscuring all text and graphics that would otherwise be present in the software specification document.

(b)(4)



Software Specification – (b)(4)

(b)(4)





(b)(4)





Software Specification – (b)(4)

(b)(4)

A large black rectangular redaction box covers the entire central portion of the page, obscuring all text and graphics that would otherwise be present in the software specification document.

(b)(4)



(b)(4)



(b)(4)



(b)(4)



Software Specification – (b)(4)

(b)(4)





(b)(4)

A large black rectangular redaction box covers the majority of the page's content, starting below the header and ending above the footer.




(b)(4)



(b)(4)




Software Specification – (b)(4)



(b)(4)

(b)(4)





(b)(4)

(b)(4)



(b)(4)



(b)(4)





(b)(4)




(b)(4)



Software Specification – (b)(4)

(b)(4)





(b)(4)



(b)(4)





Software Specification – (b)(4)

(b)(4)



(b)(4)



(b)(4)





(b)(4)







Software Specification – (b)(4)

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the section header and extending nearly to the bottom of the page.

(b)(4)



(b)(4)



## Traceability Analysis

### 1. Introduction

The development of the 3Shape Abutment Designer™ Software implements Traceability Analysis by a Trace Matrix (see *Trace Matrix* below).

**Note:** The 3Shape Abutment Designer™ Software is an add on to the 3Shape Dental System™ Software (classified as 872.3661 Optical Impression Systems for CAD/CAM, Product Code NOF, 510(K) Exempt).

The trace shown is extracted from the Dental System™ Software Trace Matrix and only relates to the Abutment Designer.

The Trace Matrix links Identified Hazards, User Needs, Requirements, Specifications, and Validation and Verification.


**Note:** 3Shape Software Product Development Documentation is structured into a suite of Requirement Specifications (RS) and a suite of Software Specifications (SS).

Each Requirement Specification and Software Specification contains an embedded test protocol for Validation and Verification respectively.

Therefore, the Validation and Verification testing requirements are implicitly identified by the unique RS and SS document IDs listed in in the Trace Matrix.

TITLE: DENTAL SOFTWARE SPECIFICATION - ARCHITECTURE

## 2. Trace Matrix



TITLE: DENTAL SOFTWARE SPECIFICATION - ARCHITECTURE

(b)(4)





TITLE: DENTAL SOFTWARE SPECIFICATION - ARCHITECTURE

(b)(4)





3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

## Software Development Environment Description

(b)(4)













































































## Human Factors Engineering Checklist

### MANUFACTURER

**Manufacturer's Name and Address:**

**3Shape A/S, Holmens Kanal 7, 1060 Copenhagen K, Denmark**

**Type and Name(s) of Medical Device(s):**

**3Shape Abutment Designer™ Software**

Element	Covered In
1 Introduction 1.1 Use related Hazards 1.2 Use Scenarios Resulting in Hazards	N/A













(b)(4) [REDACTED] 0

TITLE: RELEASE PROTOCOL - (b)(4) [REDACTED]

(b)(4) [REDACTED]

























































Validation Result – (b)(4)

(b)(4)


A large black rectangular redaction box covers the majority of the page content, starting below the section header and ending above the footer. The text "(b)(4)" is written in red at the top left corner of this redacted area.



Validation Result - (b)(4)

(b)(4)







Validation Result – (b)(4)

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the section header and extending nearly to the bottom of the page.

(b)(4)



(b)(4)





Requirement Specification - (b)(4)

(b)(4)




(b)(4)





Validation Result - (b)(4)

(b)(4)



(b)(4)



(b)(4)



(b)(4)




(b)(4)





Validation Result – (b)(4)

(b)(4)





(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)




(b)(4)



(b)(4)








(b)(4)







Validation Result - (b)(4)

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the section header and extending nearly to the bottom of the page.

(b)(4)





Requirement Specification - (b)(4)


(b)(4)






(b)(4)

A large black rectangular redaction box covers the majority of the page's content, starting below the header and ending above the footer.

Verification Result – (b)(4) 

(b)(4) 


(b)(4)






(b)(4)






(b)(4)



Verification Result – (b)(4) 

(b)(4) 

(b)(4)



(b)(4)



(b)(4)



(b)(4)







Verification Result – (b)(4) [Redacted]

(b)(4)

[Large redacted area covering the majority of the page content]

(b)(4)




(b)(4)



Verification Result – (b)(4) 

(b)(4)






(b)(4)

(b)(4)



Verification Result – (b)(4) 

(b)(4)



(b)(4)





(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)







(b)(4)



(b)(4)




Verification Result – (b)(4) 

(b)(4) 

(b)(4)







(b)(4)

(b)(4)



Verification Result – (b)(4) 

(b)(4) 

(b)(4)






(b)(4)



(b)(4)





(b)(4)

(b)(4)



(b)(4)



(b)(4)





Veirfication Result - (b)(4)

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the section header and extending nearly to the bottom of the page.

(b)(4)





(b)(4)



(b)(4)





Software Specification – (b)(4)

(b)(4)



(b)(4)



(b)(4)




(b)(4)



(b)(4)










Verification Result—(b)(4) R



(b)(4)



(b)(4)





3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

## Revision Level History

(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)







3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

## Unresolved anomalies

The list of the 3 remaining software anomalies in 3Shape Abutment Designer™ Software (b)(4) can be seen in

Table 1.

(b)(4)

**Table 1: List of open bugs in 3Shape Abutment Designer™ Software 2015-1 (version 2.15.2.0)**

Bug ID	Description	Safety Impact	Resolution	Justification for Resolution
--------	-------------	---------------	------------	------------------------------

(b)(4)

(b)(4)





3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

---

## Labeling

3Shape Abutment Designer™ is an add-on module to 3Shape Dental System™ which is classified as 872.3661 Optical Impression Systems for CAD/CAM, Product Code NOF (510(K) Exempt). Throughout this submission there will be references to 3Shape Dental System™, but only the documentation relevant to 3Shape Abutment Designer™ has been included.

This VOL\_013 includes the following documents which are excerpts of the 3Shape Dental System™ labeling:

3Shape Dental System Software Brochure  
Manual DS-2.15.2.0-A-En version

appendix 13-1  
appendix 13-2

3shape 

# Dental System™

Industry-leading scanning and CAD solutions



Excerpt  
version

## Get better esthetics and save costs with customized abutments

More and more labs are customizing abutments because they get improved clinical results, better esthetics and cost savings. Abutment Designer™ lets you automatically design the customized abutment and emergence profile with smooth transitions and optimal esthetics.

**COMING Auto-Abutment**  
Automatically generate an abutment that optimally fits the crown design.

**All types of abutments**  
Create customized abutments, screw-retained crowns and anatomical abutments.

**Easily adjust to Gingiva**  
Snap emergence profile to gingiva or use freeform tools.

**COMING Bone information**  
Improve abutment design by visualizing the jawbone.

**Abutment bevel**  
Easily add abutment bevel to your design for optimal crown support.

**NEW Add abutment position guides**  
Help your dentist easily and accurately place abutments, temporaries or veneers using custom designed positioning guides created in the same workflow.


AVINENT™  
DENSPLY BIOMET 3i™  
DIO® Nobel Biocare phibo® camlog  
straumann ZIMMER | dental  
Neo Biotech BIOHORIZON


**More than 50 implant libraries**  
Get abutments from both original and compatible-with implant manufacturers covering all global and key regional players.  
[www.3shapedental.com/implant-systems\\*](http://www.3shapedental.com/implant-systems*)

**COMING Bone information for better clinical results**  
Include bone data from CBCT scans to create better abutments without collisions and optimally shaped to the thickness of the gingiva.

## When accuracy and esthetics matter


Expand your lab's business by offering a high value service that is rapidly growing in demand. Design sophisticated implant bars and bridges with the utmost precision using flexible tools and ISO-documented 3Shape scanner accuracy.





**Design any type of bar**  
Support for a wide range of standard bar and bridge types such as Primary, Dolder, Hader, Hybrid, Canada, and Wrap-around.

**Extensive range of implant libraries\***  
Support for direct-seating and multi-unit abutments.




**Proven bar accuracy with 3Shape scanners**  
Accuracy on all 3Shape scanners is verified and documented according to ISO12836 and specially designed implant objects.

**NEW Freely apply different types of bars**  
Combine multiple types in a single design with the new smooth surface transition that ensures an optimal finish.


**NEW Improved finish with new Cut-to-gingiva**  
Match the bar to the underlying gingiva and achieve high quality with the newly enhanced surface finish.

**Add any type of attachment**  
Freely add slide or ball attachments, locators and retention holes. Adjust position and angle if needed.

\*Only libraries using FDA-cleared components permitted in the US.



**COMING Perfectly matching secondary structures**  
Design your bar and then easily and precisely create its matching secondary according to the bar's attachments and form.




**Esthetic implant bridges with gingiva**  
Easily design implant bridges complete with gingiva based on the final anatomy and a gingiva boundary that you define.



**COMING Designs obtained directly from the Denture or Wax-up**  
Scan the denture and virtually cut back to facilitate design of an optimal implant bridge, or scan an original wax-up bridge to create a file for copy milling.

## 3Shape develops 3D technologies for dental practices and labs



---

## Backing labs with care, technology and expertise



3Shape LABcare™ is an integral part of your annual 3Shape subscription. It bundles services that are designed to ensure your investment, secure maximum uptime, and help you get the most from your solution year after year.

Find 3Shape online



Questions? Contact FDA/CDRH/OCE/DID at [CDRH-FOISTATUS@fda.hhs.gov](mailto:CDRH-FOISTATUS@fda.hhs.gov) or 301.796.8118

# 3Shape Dental System 2015

## User Manual

(b)(4) Draft User Manual





























































































































## Appendix A: System Requirements

Item	Minimum Requirements	Recommended
<b>OS</b>	Windows 7 32-bit Professional*	Windows 7 64-bit Professional
<b>RAM</b>	4GB	8GB (16GB**)
<b>Video Card</b>	512MB DirectX 10 (1GB DirectX 10*) NVIDIA GeForce	1GB DirectX 11 (2GB DirectX 11*) NVIDIA GeForce
<b>Available HDD Space</b>	250GB	500GB (1TB***)
<b>CPU</b>	Intel Core i5 or equivalent	Intel Core i7 or equivalent
<b>Monitor resolution</b>	1920 x 1080 pixels	1920 x 1080 pixels
<b>3D Mouse</b>		3DConnexion SpaceMouse™ Pro
<b>Network</b>	Internet connection	
<b>USB ports</b>	USB 2.0 port for 3Shape desktop scanner****	
<b>Mouse</b>	Mouse with wheel button support	
<p>* Minimum requirement for the Implant Studio is Windows 7 64-bit with 1GB DirectX 10 video card. 2GB DirectX 11 video card is recommended.</p> <p>** For simultaneous scanning and modelling of large cases, we recommend 16GB RAM.</p> <p>*** We recommend 1TB Hard Drive if used as a stand-alone system or a server with the order folder.</p> <p>**** D250 scanners require serial port. For D250 scanners, 3Shape supports only PCs purchased from 3Shape as very few PCs are compatible.</p>		

### Note!



- 3Shape Desktop scanners work only on USB 2.0 ports.
- 3Shape Desktop Scanners do not work on shared USB connections.
- It is recommended to connect keyboard, mouse and dongle to free USB 3.0 ports, when available.



**Note!** For D250 scanners, we recommend keeping the existing PC delivered with the system for scanning, and installing a new PC (Windows 7, 64-bit, 6-8 GB RAM) as a design station.

## Appendix G: Implant Libraries for 3Shape Dental System™

Providers of Implant Solutions	Implant Systems							Implant Solutions				Manufacturing	Consumables (for local use)				Additional info			
	Nobel Biocare	Biomet3i	Straumann	Zimmer Dental	Astra Tech DENTSPLY	Friadent DENTSPLY	CAMLOG	Other systems	1-Piece (Titanium)	1-Piece (Zirconia)	2-Piece (Ti base + Zr-Abutm)		Implant Bars & Bridges	Digital implant models	Providers milling center	Local milling by 3rd Party		Blanks with pre-milled interface	Titanium bases	Intraoral Scan bodies
<b>Alpha Bio-Tec. (IL)</b> www.alpha-bio.net							- Alpha Bio Tec	X	X	X	X		X	X	X					
<b>Argen (US)</b> www.argen.com	+	+	+	+	+	+		X	X			X				X	X	X	X	- FDA, CE
<b>ATLANTIS™ (SE)</b> www.dentsplyimplants.com	+	+	+	+	+	+	- Keystone Dental - BioHorizons	X	X			X								- FDA, CE
<b>BEGO (DE)</b> www.bego.com	+	+	+	+	+	+	- BEGO Implant Systems - and more	X	X	X	X	X					X	X		- FDA, CE - CE-labelled prosthetic screws
<b>BioComp (NL)</b> www.biocomp.eu							- BioComp	X	X				X	X	X	X	X	X		- CE, ISO
<b>Biodenta (CH)</b> www.biodenta.com	+	+	+	+	+	+	- Biodenta	X	X	X	X	X						X	X	- FDA, CE pending
<b>BioHorizons (US)</b> www.biohorizons.com	+		+	+			- BioHorizon	X	X				X	X	X	X				- FDA, CE
<b>Biomet (US)</b> www.biomet.com	+	+	+					X	X	X		X						X		- FDA, CE - Encode® healing abutments
<b>Biotech International (FR)</b> www.biotech-dental.com							- Biotech			X	X		X			X	X	X		- CE
<b>Bredent (DE)</b> www.bredent.com							- Sky			X			X			X	X	X		- FDA, CE





<b>CADBLU (US)</b> www.cadbludental.com	+	+	+	+	+	- Bio Horizons - MIS Implants	X		X	X	X	X					X	X	- FDA
<b>CAMLOG (CH)</b> www.camlog.com						+ - CONELOG - iSy	X		X	X		X	X		X	X			
<b>CAP (US)</b> www.cap-us.com	+	+	+	+	+	+		X		X	X	X	X	X	X	X		X	
<b>CMC (US)</b> www.custom-milling.com	+	+	+	+	+	+	+	- Bio Horizons - Osstem - LifeCore - Ultra-Lock - Neoss	X	X	X	X		X				X	- FDA, CE
<b>Core3D International (NL)</b> www.core3dcenters.com	+	+	+	+	+	+	+	- Avinent - BioComp - BioHorizons - MIS Implants - and many more	X		X	X	X	X	X	X	X	X	- FDA, CE, Health Canada, TGA (AUS), Taiwan regulatory - MHLW: Ministry of Health, Labour and Welfare (JP)
<b>C-Tech Implant (IT)</b> www.c-tech-implant.com							- C-Tech Implant	X		X	X			X		X		- CE	
<b>Degudent (DE)</b> www.degudent.com	+	+	+	+	+	+	+	- Medentika	X	X	X			X	X	X			
<b>Dental Consulting (DE)</b> www.gadau-consulting.com	+		+	+	+	+	+	- ICX-Medentis - Osstem	X		X	X	X	X	X	X	X	X	- FDA pending, CE
<b>Dentaurum Implants (DE)</b> www.dentaurum-implants.de								- Dentaurum Implants			X				X		X	X	- CE, PAL
<b>Dentegris Deutschland (DE)</b> www.dentegris.de								- Dentegris Implants			X	X			X		X		
<b>Dentsply-Friident (DE)</b> www.dentsply-friident.com							+	- Xive - Ankylos	X	X	X			X				- CE	
<b>DESS (ES)</b> www.dess-abutments.com	+	+	+	+	+	+		- BioHorizon	X		X			X	X		X	X	- CE - ISO 13485 - ISO 9001
<b>Digital Dental Group-DDG (ISR)</b> www.ddg-scanlab.com	+	+	+	+	+	+	+	- Microdent - BioHorizons - MIS-Implants - and more	X	X	X	X			X	X	X	X	- FDA pending, CE
<b>DIO-Implants (KR)</b> www.dioimplant.com	+							- DIO Implants	X	X	X	X		X	X			- FDA, CE, KFDA, SFDA	

<b>Elos Medtech Pinol (DK)</b> www.elosmedtech.com	+	+	+	+	+	+	+	- Neoss	X	X	X	X	X	X	X	X	X	X	- FDA, CE	
<b>EuroTeknika (FR)</b> www.euroteknika.com	+	+	+	+	+			- Anthogyr - Biotech - Easy Implant - MIS Implants - and more	X	X	X	X	X		X	X	X	X	- FDA, CE	
<b>GC Advanced Technologies (US)</b> www.gc-at.com	+	+	+	+	+			- BioHorizons - Sybron	X	X	X		X							
<b>Glidewell (US)</b> www.glidewell dental.com	+	+	+	+	+	+	+	- Prisma - DentalCraft - Keystone Dental - Neoss	X	X	X	X	X			X	X	X	- FDA, CE pending	
<b>Heraeus Kulzer (DE)</b> www.heraeus.com	+	+	+	+	+	+	+	- Thommen	X	X	X		X	X						
<b>Ivoclar Wieland (DE)</b>	+	+	+	+	+	+		- Denta - Aurum Implants	X		X	X		X	X		X	X		
<b>LaStruttura (IT)</b> www.lastruttura.it	+	+	+	+	+	+		- Sweden&Martin - Megagen - Prodent - and many more	X	X		X	X	X			X	X	- CE 93/42	
<b>Medentika (DE)</b> www.medentika.de	+	+	+	+	+	+	+	- Medentika - M-Implant	X		X	X		X	X	X	X	X	- CE	
<b>Medentis Medical (DE)</b> www.medentis.de	+		+		+	+		- ICX - Templat	X		X	X	X	X			X	X	- CE	
<b>Medical Production (FR)</b> www.medical-production.eu						+		- Euroteknika	X		X	X	X		X	X	X	X	- CEO 499 - ISO 9001 - ISO 13485	
<b>MIS Implants Technology (IL)</b> www.mis-implants.com								- MIS Implants			X				X		X		- FDA, CE	
<b>Neodent (BR)</b> www.neodent.com.br	+		+					- Neodent - Implant Systems	X	X	X	X			X			X		
<b>Neoss (UK)</b> www.neoss.com	+	+	+	+	+			- Neoss	X		X	X	X	X	X	X	X		X	- FDA, CE
<b>Nobel Biocare (CH)</b> www.nobelbiocare.com	+	+	+	+	+				X	X	X			X					- FDA, CE	
<b>NT Trading (DE)</b> www.nt-trading.com	+	+	+	+	+	+	+	- Thommen - Sweden&Martin - a	X		X	X		X	X	X	X	X	X	- FDA, CE, - CMDCAS, GOST - R, TGC - Mexico

																	- BEGO									Certification			
<b>Phibo (ES)</b> www.phibo.com	+	+	+	+	+	+	+	+	+								- PHIBO - BTI - Sweden&Martina - MIS Implants - and many more	X	X	X	X	X	X				X	X	- FDA, CE - ISO 13485 - ISO 9001
<b>Prowital (DE)</b> www.prowital.de																	- Prowital			X						X	X	X	- CE - EN-ISO 9001-V472008 - EN ISO-13485-2003AC2009 - Certificate Directive-93-42 EWG
<b>Ritter Implants (DE)</b> www.ritterimplants.com			+	+													- Alpha-BioTec. - MIS Implants	X		X	X			X	X	X	X	- FDA, CE	
<b>Straumann (CH)</b> www.straumann.com			+																X	X	X			X	X			- CE	
<b>Sweden&amp;Martina S.p.A (IT)</b> www.swedenmartina.com																	- Sweden&Martina	X	X	X	X	X	X			X	X	X	
<b>Target3D (USA)</b> www.target3d.com	+	+	+	+	+	+	+	+	+								- Osstem - Bio Horizons - MIS-Implants - and more	X	X	X	X	X		X	X	X	X	X	- FDA, CE
<b>Thommen Medical (CH)</b> www.thommenmedical.com																	- Thommen			X	X					X	X	- FDA, CE	
<b>TRI Implants (CH)</b> www.tri-implants.com																	- TRI Implants			X						X	X	X	
<b>3dental (ES)</b> www.3dental.es	+	+	+	+	+	+											- BioComp - Sweden&Martina - DYNA	X	X		X	X	X				X		

+ -Available for Dental System™      x -Available from the manufacturer

 **Note!** All information is given without guarantee and based exclusively on information made available by the implant system providers. Please contact your local 3Shape re-seller to obtain the most current list.

 **Note!** Only Implant Libraries based on Dental Implants with 510(k) clearance can be used in Abutment Designer™ in the United States. Cleared libraries must be activated by 3Shape. Please ask your Implant Library Provider to contact 3Shape if you wish to use a cleared library that has not already been activated. The software will block any attempt to use a non-cleared library.

## Appendix H: Terms and Abbreviations

The following table explains file type abbreviations used in 3Shape Dental System.

File Format	Description
3ML	Zipped, compressed XML file typically used for setup and customization. 3Shape proprietary format.
3SE	3Shape and Sirona proprietary format for exported orders. Orders from Sirona system, can be imported into 3Shape system.
3OX	3Shape Communicate Order Exchange format used by partners and customers to retrieve and send orders. 3Shape proprietary format.
3OXZ	Zipped archive containing 3OX file and references to DCM models. 3Shape proprietary format.
3SI	3Shape and Sirona proprietary format for imported orders. Orders exported from 3Shape system, can be imported into Sirona system.
DCM	Dental Compressed Model file. Contains compressed 3D model data, attached objects (splines, annotations, etc.), marks and additional string properties. 3Shape proprietary format used for scans and CAD designs.
DLL	Dynamic-link library, Microsoft shared library concept.
DME	Dental System Material Export file. Contains materials, references to materials and external files. The file can be imported into another 3Shape Dental System. 3Shape proprietary format.
STL	Describes surface geometry of three-dimensional objects. Used for scans and CAD designs. Industry standard.
ULDC	3M Lava proprietary format. Order files from 3M scanners can be imported into 3Shape System.
XML	Extensible Markup Language, used for configuration files, etc.

## Appendix I : Contact Information

### **3Shape Headquarters**

Europe, Middle East &  
Africa Sales  
Holmens Kanal 7  
1060 Copenhagen K  
Denmark

**P:** +45 70 27 26 20

### **3Shape North America**

North American Sales  
Somerset Hills  
Corporate Center  
10 Independence  
Boulevard,  
Suite 150  
Warren, New Jersey  
07059, USA

**P:** +1 908 867 0144

### **3Shape (Shanghai) Co., Ltd**

Asian Sales  
Room 906, Tower A of  
Eton Place  
No. 69, Dongfang  
Road  
200120 Shanghai,  
China

**P:** +86 21 5835 2281

### **3Shape South America**

Latin American and  
Caribbean Sales  
Carrera 13 # 82-91  
Oficina 401  
110221 Bogotá,  
Colombia

**P:** +57 1 691 95 08



K151455/S001  
3shape

Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center - WO66-G609  
10903 New Hampshire Avenue  
Silver Spring, Maryland 20993-0002  
USA

FDA/CDRH/DCC

JAN 11 2016

RECEIVED

January 8, 2016

Re: 510(k) # K151455/S001 Additional Information (AI) Request -  
Replacement eCopy

Att.: Document Mail Clerk

Enclosed please find our reply to the Additional Information request received by email July 31, 2015. The following documentation has been included:

Document ID	Title	Starting Page Number
001	AI1 Reply - Cover Letter	1
002	AI1 Reply	3
003	510(k) Summary	19
004	Substantial Equivalence Comparison	23
005	Software Description	26
006	Vertical Offset Description	34
007	Device Hazard Analysis	35
008	Cyber Security Analysis	36
009	Architecture Design Chart	37

14

PAGE 1

Document ID	Title	Starting Page Number
010	Software Design Specification	43
011	Traceability Analysis	44
012	Unresolved Anomalies	49
013	Proposed Labeling	51
014	Indications for Use Statement	52
015	SW Requirement Specifications	53
016	Verification and Validation	54
017	Product Risk Assessment	57
018	SW-SOP-0011 – Handling of FDA Cleared Implant Libraries in 3Shape Software	60
019	Appendix 1 – Library Confirmation Letters	61
020	Appendix 2 - User Manual / Labeling	63

**Submission copies:** 2 (eCopy and paper). The eCopy is an exact duplicate of the paper copy.

Best regards,  
3Shape A/S



Hanne Nielsen  
Regulatory Affairs Manager

*K151455/S001*

FDA CDRH DMC

**3shape** 

DEC 09 2015

Received

Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center - WO66-G609  
10903 New Hampshire Avenue  
Silver Spring, Maryland 20993-0002  
USA

December 7, 2015

Re: 510(k) # K151455/S001


Att.: Document Mail Clerk

Enclosed please find revised electronic copy of Appendix 1 of the above mentioned submission. No changes have been made to the paper version of the documentation submitted October 14<sup>th</sup>. The following documentation has been included:

Document ID	Title	Starting Page Number
001	AI1 Reply - Cover Letter	1
019	Appendix 1 - Library Confirmation Letters	61

**Submission copies:** 2 (eCopy and paper). The eCopy is an exact duplicate of the paper copy.

Best regards,  
3Shape A/S



Hanne Nielsen  
Regulatory Affairs Manager

54  
PAGE 1



Received

OCT 19 2015

3shape 

FDA CDRH DMC

Food and Drug Administration  
 Center for Devices and Radiological Health  
 Document Mail Center - WO66-G609  
 10903 New Hampshire Avenue  
 Silver Spring, Maryland 20993-0002  
 USA

FDA CDRH DMC

OCT 19 2015

Received

October 15, 2015

*K151455/5001*

Re: 510(k) # K151455 Additional Information (AI) Request

Att.: Document Mail Clerk

Enclosed please find our reply to the Additional Information request received by email July 31, 2015. The following documentation has been included:

Document ID	Title	Starting Page Number
001	AI1 Reply - Cover Letter	1
002	AI1 Reply	3
003	510(k) Summary	19
004	Substantial Equivalence Comparison	23
005	Software Description	26
006	Vertical Offset Description	34
007	Device Hazard Analysis	35
008	Cyber Security Analysis	36
009	Architecture Design Chart	37
010	Software Design Specification	43

*68*



Document ID	Title	Starting Page Number
011	Traceability Analysis	44
012	Unresolved Anomalies	49
013	Proposed Labeling	51
014	Indications for Use Statement	52
015	SW Requirement Specifications	53
016	Verification and Validation	54
017	Product Risk Assessment	57
018	SW-SOP-0011 - Handling of FDA Cleared Implant Libraries in 3Shape Software	60
019	Appendix 1 - Library Confirmation Letters	61
020	Appendix 2 - User Manual / Labeling	63

**Submission copies:** 2 (eCopy and paper). The eCopy is an exact duplicate of the paper copy.

Best regards,  
3Shape A/S

A handwritten signature in black ink, appearing to read "Hanne Nielsen".

Hanne Nielsen  
Regulatory Affairs Manager

Food and Drug Administration  
 Center for Devices and Radiological Health  
 Document Mail Center - WO66-G609  
 10903 New Hampshire Avenue  
 Silver Spring, Maryland 20993-0002  
 USA

January 8, 2016

Re: 510(k) # K151455/S001 Additional Information (AI) Request –  
 Replacement eCopy

Att.: Document Mail Clerk


Enclosed please find our reply to the Additional Information request received by email July 31, 2015. The following documentation has been included:

Document ID	Title	Starting Page Number
001	AI1 Reply – Cover Letter	1
002	AI1 Reply	3
003	510(k) Summary	19
004	Substantial Equivalence Comparison	23
005	Software Description	26
006	Vertical Offset Description	34
007	Device Hazard Analysis	35
008	Cyber Security Analysis	36
009	Architecture Design Chart	37

Document ID	Title	Starting Page Number
010	Software Design Specification	43
011	Traceability Analysis	44
012	Unresolved Anomalies	49
013	Proposed Labeling	51
014	Indications for Use Statement	52
015	SW Requirement Specifications	53
016	Verification and Validation	54
017	Product Risk Assessment	57
018	SW-SOP-0011 - Handling of FDA Cleared Implant Libraries in 3Shape Software	60
019	Appendix 1 - Library Confirmation Letters	61
020	Appendix 2 - User Manual / Labeling	63

**Submission copies:** 2 (eCopy and paper). The eCopy is an exact duplicate of the paper copy.

Best regards,  
3Shape A/S



Hanne Nielsen  
Regulatory Affairs Manager



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455.S001 ACK Letter

**510(k) # K151455**

Reply to Additional Information #1 request dated July 31, 2015

(b)(4)





3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455.S001 ACK Letter

(b)(4)

A large, solid black rectangular redaction box covers the majority of the page's content, starting below the header and ending above the footer. The text "(b)(4)" is printed in red at the top left corner of this redacted area.



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455.S001 ACK Letter

(b)(4)

A large, solid black rectangular redaction box covers the majority of the page's content, starting below the header and ending above the footer. The text "(b)(4)" is printed in red at the top left corner of this redacted area.



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455.S001 ACK Letter

(b)(4)

A large, solid black rectangular redaction box covers the majority of the page's content, starting below the header and ending above the footer. The text "(b)(4)" is printed in red at the top left corner of this redacted area.





3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455.S001 ACK Letter

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the header and ending above the footer. The text "(b)(4)" is written in red at the top left corner of this redacted area.



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455.S001 ACK Letter

(b)(4)





3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455.S001 ACK Letter

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the header and ending above the footer. The text "(b)(4)" is written in red at the top left corner of this redacted area.

CONFIDENTIAL

Page 9



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455.S001 ACK Letter

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the header and ending above the footer. The text "(b)(4)" is written in red at the top left corner of this redacted area.



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455.S001 ACK Letter

(b)(4)

A large, solid black rectangular redaction box covers the majority of the page's content, starting below the header and ending above the footer. The text "(b)(4)" is printed in red at the top left corner of this redacted area.



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455.S001 ACK Letter

(b)(4)





3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455.S001 ACK Letter

(b)(4)





3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455.S001 ACK Letter

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the header and ending above the footer. The text "(b)(4)" is written in red at the top left corner of this redacted area.





3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455.S001 ACK Letter

(b)(4)

A large, solid black rectangular redaction box covers the majority of the page content, starting below the header and ending above the footer. The text "(b)(4)" is written in red at the top left corner of this redacted area.



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455.S001 ACK Letter

(b)(4)

A large, solid black rectangular redaction box covers the majority of the page's content, starting below the header and ending above the footer. The text "(b)(4)" is written in red at the top left corner of this redacted area.



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455.S001 ACK Letter

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the header and ending above the footer.

CONFIDENTIAL

Page 17



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455.S001 ACK Letter

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the header and ending above the footer. The text "(b)(4)" is written in red at the top left corner of this redacted area.

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

**510(K) SUMMARY – Traditional 510(K)**

**Submitter Information**

A Company Name: 3Shape A/S  
B Company Address: Holmens Kanal 7  
DK-1060 Copenhagen K  
C Company Phone: +45 7027 2620  
Company Fax: +45 7027 2621  
D Contact Person: Hanne Nielsen  
Regulatory Affairs Manager  
E Date Summary Prepared: October 12, 2015

**Device Identification**

A Trade/proprietary Name: 3Shape Abutment Designer™  
Software  
B Common Name: Abutment Designer  
C Device Classification Name: Endosseous Dental Implant Abutment  
C Regulation Number: 872.3630  
C Classification: Class II  
D Product Code: NHA

**Predicate Device**

The 3Shape Abutment Designer™ Software has equivalent intended use and technical characteristics as the Sirona Dental CAD/CAM System (K100152).

3Shape Abutment Designer™ Software does not provide any physical parts that can come into contact with the patient.

The software output is restricted in the U.S. to be manufactured by a holder of an implant abutment 510(k) or milled per the specific instructions provided by the holder of a 510(k) marketing clearance for a patient specific implant abutment.

### 3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

This is indicated in the device labelling.

Therefore, the differences between the Device and the predicates do not raise additional concerns concerning the Device's safety and effectiveness.

#### **Intended Use**

The 3Shape Abutment Designer™ Software is intended for designing endosseous dental implant abutments based on digital scans of 3D models.

#### **Device Description**

The 3Shape Abutment Designer™ Software receives scan data containing topographical characteristics of real teeth, dental impressions, or stone models. The software provides the user with the ability to create matching endosseous dental implant abutments using computer aided design.

The output of the device is a computer file containing the abutment(s) in digital form which can be used by manufacturers that hold an implant abutment 510(k) or are milling per the specific instructions provided by the holder of a 510(k) marketing clearance for a patient specific implant abutment.

The 3Shape Abutment Designer™ Software can be run on properly configured "off-the-shelf" PC hardware running Microsoft® Windows and using a standard keyboard and mouse.

#### **Scientific Concept**

The underlying scientific concept is the use of CAD/CAM technology for the optical acquisition of the topographical characteristics of dental impressions, and models and the design of individual mesostructures using recorded data (CAD).

#### **Summary of the technological characteristics**

The 3Shape Abutment Designer™ is a software only device programmed in the Delphi programming language and has the following PC/laptop requirements:

## 3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

Item	Minimum Requirements	Recommended	Predicate Device
<b>OS</b>	Windows 7 32-bit Professional*	Windows 7 64-bit Professional	Windows Vista 32-bit
<b>RAM</b>	4GB	8GB (16GB)	6GB
<b>Video Card</b>	512MB DirectX 10 (1GB DirectX 10) NVIDIA GeForce	1GB DirectX 11 (2GB DirectX 11) NVIDIA GeForce	512 MB DirectX 10 NVIDIA Quadro
<b>Available HDD Space</b>	250GB	500GB (1TB)	500GB
<b>CPU</b>	Intel Core i5 or equivalent	Intel Core i7 or equivalent	Intel Core i7
<b>Monitor resolution</b>	1440 x 900 pixels	1920 x 1080 pixels	Unknown
<b>3D Mouse</b>	None	3DConnexion SpaceMouse™ Pro	Unknown
<b>Network</b>	Internet connection		Unknown
<b>Ports</b>	USB 2.0 port for 3Shape desktop scanner		Serial Port PCI card for inEOS scanner
<b>Mouse</b>	Mouse with wheel button support		Unknown

### Nonclinical Testing

Software, hardware, and integration verification and validation testing was performed in accordance with the FDA Guidance Document "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices" (Issued on May 11, 2005).

Prior to release, verification and validation testing of the 3Shape Abutment Designer™ Software has been completed using the approved acceptance criteria: Each user need has its own validation acceptance criteria; each specification has its own verification acceptance criteria; bug verification consists in ensuring issue is not reproducible; issues reported by beta partners must be reviewed and handled appropriately; beta partners must report functionality to be an improvement over previous version.

The validation suite includes validation of implemented mitigations related to device hazards identified in the risk management procedures.

All test results have been reviewed and approved, showing the 3Shape Abutment Designer™ Software to be safe and effective.

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

**Clinical Testing**

Clinical testing is not a requirement and has not been performed.

**Conclusion**

Based on a comparison of intended use, indications, construction materials, principle of operations, features and technical data, and the test results, the 3Shape Abutment Designer™ Software is found to be substantial equivalent with the Predicate Devices.



## Substantial Equivalence Comparison

### 1. Predicates

The 3Shape Abutment Designer™ Software has the same intended uses and technical characteristics as the Sirona Dental CAD/CAM System (K100152) as listed in “Table 1: Predicate”

**Table 1: Predicate**

<b>Predicate</b>	<b>Manufacturer</b>	<b>510(k) number</b>	<b>Product code</b>
Dental CAD/CAM System	Sirona	K100152	NHA*

\* Endosseous dental implant abutments, 21CFR872.3630

## 2. Intended Use Comparison

### 2.1. 3Shape Abutment Designer™ Software

The Device's Intended Use, Intended users, and Intended Operational is stated in the "VOL\_001\_Administrative Documents" volume of this submission and reproduced here:

*The 3Shape Abutment Designer Software is intended as an aid the restoration of chewing function in partially or fully edentulous mandibles and maxillae. The 3Shape Abutment Designer Software is intended for use by a dental practitioner or dental laboratory staff for designing the patient specific component of a two-piece, one-piece, or hybrid dental implant abutment. The single or multi-unit abutment design is intended to be used by the manufacturer of an endosseous dental implant abutment to create the final device.*

*Intended users are Dental Practitioners and Dental Laboratory staff.*

*Intended Operational Environment is Dental Laboratories."*

In the following sections, the similarities with the predicate are discussed.

### 2.2. Dental CAD/CAM System

The predicate's Intended Use can be extracted from the FDA 510(k) Premarket Notification Database:

*"The Sirona Dental CAD/CAM System is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations. The system consists of three major parts: TiBase, InCoris mesostructure, and CAD/CAM software. Specifically, the InCoris mesostructure and TiBase components make up a two-piece abutment which is used in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity.*

*The InCoris mesostructure may also be used in conjunction with the Camlog Titanium base CAD/CAM (types K2244.xxxx) (K083496) in the Canmlog Implant System. The CAD/CAM software is intended to design and fabricate the InCoris mesostructure. The InCoris mesostructure and TiBase two-piece abutment is compatible with the following implants systems:*

*[List of compatible implant systems omitted for brevity]*

Underlined segments indicate similarity with 3Shape Abutment Designer™ Software.

Note, the predicate is a CAD/CAM System bundled with physical two-piece abutments; the 3Shape Abutment Designer™ Software is CAD/CAM only.

### 3. Characteristics Comparison

A Characteristics Comparison can be seen in "Table 2: Substantial Equivalence Chart".

**Table 2: Substantial Equivalence Chart**

<b>Feature name</b>	<b>3Shape Abutment Designer™ Software</b>	<b>Sirona Dental CAD/CAM System (K100152)</b>
Graphical UI	Yes	Yes
Windows OS platform	Yes	Yes
Uses standard PC hardware	Yes	Yes
Digitally imports topography of teeth by 3D Scan	Yes	Yes
Uses 3D CAD design tools	Yes	Yes
Custom abutment design	Yes	Yes
Screw retained design	Yes	Yes
Implant Bar design	Yes	Yes
Export to remote milling machine by internet	Yes	Yes
Network Protocol	Internet/TCP-IP	Internet/TCP-IP
Intended users	Dental practitioners and dental labs	Dental practitioners and dental labs
Output type	Computer file	Computer file
Device submission includes pre-manufactured prosthetics*	No	Yes

\* Endosseous dental implant abutments as per 21CFR872.3630

### 4. Conclusion

The 3Shape Abutment Designer™ Software and the predicate only deviate significantly in the cases where the predicate is bundled with a physical dental implant abutment. 3Shape Abutment Designer™ Software does not provide any physical parts that can come into contact with the patient.

The software output is restricted in the U.S. to be manufactured by a holder of an implant abutment 510(k) or milled per the specific instructions provided by the holder of a 510(k) marketing clearance for a patient specific implant abutment.

This is indicated in the device labelling.

The differences between the 3Shape Abutment Designer™ Software and the predicate do not raise additional concerns with respect to the safety and effectiveness of the 3Shape Abutment Designer™ Software.

Based on the information presented, we conclude Substantial Equivalence between the predicate and the 3Shape Abutment Designer™ Software.



3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

---

## Software Description

(b)(4)





3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

## **2. Intended Use and Operational Environment**

The 3Shape Abutment Designer Software is intended as an aid the restoration of chewing function in partially or fully edentulous mandibles and maxillae. The 3Shape Abutment Designer Software is intended for use by a dental practitioner or dental laboratory staff for designing the patient specific component of a two-piece, one-piece, or hybrid dental implant abutment. The single or multi-unit abutment design is intended to be used by the manufacturer of an endosseous dental implant abutment to create the final device.

Intended users are Dental Practitioners and Dental Laboratory staff.

Intended Operational Environment is Dental Laboratories.

## 3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

**3. Definitions**


There are four types of Abutments supported by the 3Shape Dental System: Customized, Anatomical, Screw Retained Crown, and Wax-up.

3Shape uses the term "Abutment" to denote all four types. Digitally created abutments are classified as an Endosseous dental implant abutment.

The 3Shape Abutment Designer™ only supports modification of the Abutment itself as exemplified on *Figure 1: Customized Abutment parts*.

A digital representation of the Implant is required for the software to work, but only to ensure the Abutment will fit on the Seating Geometry / Interface of the Implant. The 3Shape Abutment Designer™ does not provide any means to design, alter, or manufacture any part of the Implant including, but not limited to, the abutment-to-implant interface.

**Figure 1: Customized Abutment parts**

**Terms used**

The 3Shape Abutment Designer™ Software supports four types of Abutments – however this only has impact on the initial guess of abutment shape as well as slight differences in the User Interface for designing the Abutment.



3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

There is no difference between the four types on abutment function, manufacturing requirements, output types, etc.

The output of any abutment designed is a computer file of the surface model.

Please find a detailed description in *Table 2: Abutment Types*

**Table 2: Abutment Types**

Type	Description
Customized Abutment	<p>The Customized Abutment comes in 3 subtypes: Custom, Robotic, and Bar Interface.</p> <p>All provide the user with the ability to free-hand modify the Top Cap and Emergence Profile as depicted on <i>Figure 1: Customized Abutment parts</i>.</p> <p>The Custom Abutment initial shape guess is based on an anatomic heuristic</p>  <p>The Robotic initial shape guess is based on a classic standard abutment.</p>  <p>The Bar Interface initial shape guess is based on a cylinder.</p>

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

	
<p>Anatomical Abutment</p>	<p>This type of abutment comes with an initial shape guess that matches a crown.</p> 
<p>Screw Retained Crown</p>	<p>This type of abutment comes with an initial shape guess that can be used as a crown and has a screw-hole in the centre.</p> 
<p>Wax-up abutment</p>	<p>This type of abutment comes with an initial shape guess based on a scan of a wax-model. Otherwise, it is equivalent to Anatomical abutment.</p>

Also, throughout the submission the following terms are applicable.

**Table 1: Terms**

Term	Definition
Provider	Any creator of an Implant Library.
Implant Library	A digital representation of an Implant System
Implant System	A physical Dental Implant compatible with one of more abutments
Abutment	Endosseous dental implant abutment as regulated under 872.3630





3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

<p>Milling Center</p>	<p>An industry term for a large scale provider of Milling. I.e. manufacturer that hold an implant abutment 510(k) or are milling per the specific instructions provided by the holder of a 510(k) marketing clearance for a patient specific implant abutment.</p> <p>Typically, the actual manufacturing is geographically centralized, which makes digital transfer of CAD designed surface files convenient.</p> <p>3Shape does not own or operate any milling centers.</p>
<p>Local Milling</p>	<p>When the Dental Lab and milling machinery is co-located. 510(k) restrictions also apply as described under "Milling Center".</p> <p>3Shape does not manufacture or market milling machinery.</p>



3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

#### 4. Ensuring Regulatory Compliance of Abutments

Abutments are designed to interface with Dental Implants and 3Shape Abutment Designer™ requires a digital representation of the Dental Implant in order to create a functional Abutment.

The digital representations of Dental Implants (called Implant Systems) is collected in files called Implant Libraries.

In the U.S. only Abutments designed against Dental Implants with a 510(k) clearance are allowed in the 3Shape Abutment Designer™ (see *Design stops in the software* below). All settings on the Implant systems will be locked for editing.

Providers of Implant Libraries to be used in the U.S. must supply 3Shape with written documentation in order to have their Libraries activated in the software.

It is the responsibility of the Provider to ensure that the digital representation correctly represent the cleared physical parts.

##### Design stops in the software

The software provides a “wizard”-type functionality that enable Implant System Providers to create Implant Libraries to be used in the 3Shape Abutment Designer™ Software.

The Provider can subsequently export the library to a computer file for distribution.

It is the Provider’s responsibility to ensure to lock the Implant Library. This will prevent anyone but the Provider to modify it. In the U.S., the Provider is the 510(k) holder of the Implant System.


Additionally, Providers of Implant Libraries to be used in the U.S. must supply 3Shape with written documentation in order to have their Libraries activated in the software.

Design stops in the software are implemented according to the following four steps:

- i. Libraries are locked to prevent anyone but the Provider to modify it.
- ii. The software prevents the end-user from modifying the Implant Library after it has been imported to his or her system.
- iii. The software will block any attempt to use a non-cleared Implant System in the U.S.
- iv. The software provides design stops for: Abutment Gingival Margin Diameter, Abutment Gingival Margin Height, Abutment Total Height, Abutment Angulation, and Abutment Blank Limitation.

The following table shows the trace of Requirements, Specifications, Validation, and Verification tests relevant.

Items	Requirements	Specifications	Validation	Verification
-------	--------------	----------------	------------	--------------





3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

(b)(4)

A large black rectangular redaction box covering the majority of the page content.

## 5. Manufacturing of Abutments


(b)(4)

A large black rectangular redaction box covering the majority of the page content.

## Vertical Offset Description

Vertical offset does not influence on the distance between abutment and implant.

Vertical offset is length of additional vertical cylinder surface before emergence profile to modify the emergence exit profile.



---

## Device Hazard Analysis

(b)(4) Testing, Performance, Technical Data





3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

---

## Cyber Security Analysis

### 1. Introduction

Please refer to the Cyber Security Analysis enclosed in this volume.

### 2. Table of Contents

(b)(4)	Cyber Security Analysis
--------	-------------------------



## Cybersecurity analysis report –

TITLE (b)(4) Cybersecurity Analysis Report

(b)(4)



(b)(4)





(b)(4)



(b)(4)



FORM-0002

CONFIDENTIAL

PAGE 4 OF 10

(b)(4)

FORM-0002

CONFIDENTIAL

PAGE 5 OF 10

(b)(4)

FORM-0002

CONFIDENTIAL

PAGE 6 OF 10

(b)(4)



FORM-0002

CONFIDENTIAL

PAGE 7 OF 10

(b)(4)

FORM-0002

CONFIDENTIAL

PAGE 8 OF 10

(b)(4)



FORM-0002

CONFIDENTIAL

PAGE 9 OF 10


(b)(4)





---

## Architecture Design Chart



TITLE: ARCHITECTURE DESIGN CHART

(b)(4)



TITLE: ARCHITECTURE DESIGN CHART

(b)(4)




TITLE: ARCHITECTURE DESIGN CHART

(b)(4)



TITLE: ARCHITECTURE DESIGN CHART



TITLE: ARCHITECTURE DESIGN CHART

(b)(4)



---

## Software Design Specification

### 1. Introduction

Please refer to the Software Design Specification enclosed in this volume.

Please note that the Verification Protocol templates are included in the Software Specifications.

For links between the Requirements Specifications and the Software Specifications, please see "VOL\_008\_Traceability Analysis" enclosed in this submission.


### 2. Software Design Specification, Table of Contents

(b)(4)



Software Specification – (b)(4)

(b)(4)






























Software Specification – (b)(4)

(b)(4)

A large black rectangular redaction box covers the entire central portion of the page, obscuring all text and graphics that would otherwise be present in the software specification document.



(b)(4)

(b)(4)



(b)(4)



(b)(4)




(b)(4)







Software Specification – (b)(4)



(b)(4)





Software Specification – (b)(4)

(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)





(b)(4)





Software Specification – (b)(4)

(b)(4)



(b)(4)



(b)(4)



(b)(4)




(b)(4)



Software Specification - (b)(4)

(b)(4)



(b)(4)





(b)(4)



(b)(4)





Software Specification – (b)(4)

(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)







Software Specification – (b)(4)

(b)(4)

A large black rectangular redaction box covers the entire central portion of the page, obscuring all text and graphics that would otherwise be present in the software specification document.

(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)





(b)(4)



(b)(4)





Software Specification – (b)(4)

(b)(4)




(b)(4)




(b)(4)





(b)(4)



(b)(4)

(b)(4)





(b)(4)



(b)(4)



(b)(4)




(b)(4)



(b)(4)





Software Specification – (b)(4)

(b)(4)

A large black rectangular redaction box covers the entire central portion of the page, obscuring all text and graphics that would otherwise be present in the software specification document.

(b)(4)



(b)(4)





(b)(4)





Software Specification – (b)(4)

(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)





## Traceability Analysis

### 1. Introduction

The development of the 3Shape Abutment Designer™ Software implements Traceability Analysis by a Trace Matrix (see *Trace Matrix* below).

**Note:** The 3Shape Abutment Designer™ Software is an add-on to the 3Shape Dental System™ Software (classified as 872.3661 Optical Impression Systems for CAD/CAM, Product Code NOF, 510(K) Exempt).

The trace shown is extracted from the Dental System™ Software Trace Matrix and only relates to the Abutment Designer.

The Trace Matrix links Identified Hazards, User Needs, Requirements, Specifications, and Validation and Verification.

**Note:** 3Shape Software Product Development Documentation is structured into a suite of Requirement Specifications (RS) and a suite of Software Specifications (SS).

Each Requirement Specification and Software Specification contains an embedded test protocol for Validation and Verification respectively.

Therefore, the Validation and Verification testing requirements are implicitly identified by the unique RS and SS document IDs listed in in the Trace Matrix.

TITLE: DENTAL SOFTWARE SPECIFICATION - ARCHITECTURE

## 2. Trace Matrix

(b)(4)



TITLE: DENTAL SOFTWARE SPECIFICATION - ARCHITECTURE

(b)(4)



TITLE: DENTAL SOFTWARE SPECIFICATION - ARCHITECTURE

(b)(4)



TITLE: DENTAL SOFTWARE SPECIFICATION - ARCHITECTURE

(b)(4)





3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

## Unresolved anomalies

The list of the 3 remaining software anomalies in 3Shape Abutment Designer™ Software (b)(4) can be seen in


Table 1.

(b)(4)

Table 1: List of open bugs in 3Shape Abutment Designer™ Software (b)(4)

Bug ID	Description	Safety Impact	Resolution	Justification for Resolution
--------	-------------	---------------	------------	------------------------------

(b)(4)



## 3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

## Proposed Labelling

The 3Shape Abutment Designer™ Software is a prescription device and is exempt from needing adequate directions for lay use.

The following labelling is included (please see "VOL\_013 Labelling" of this submission):

- PDF copy of User Manual DS-2.15.2.0-A-EN
- Dental System™ Brochure

The application is for download only, and hence there is no CD label.

### Device Specific Requirements

There are no applicable requirements in a device-specific regulation.

### Special Controls Document

The 3Shape Abutment Designer™ Software falls into the following product group:

#### NHA

The guidance document entitled "Class II Special Controls Guidance Document: Root-Form Endosseous Dental Implants and Endosseous Dental Implant Abutments" serve as the special control.

The Guidance Document lists the below labelling requirements.

Requirement	Referenced in
Labelling in sufficient detail to satisfy the requirements of 21 CFR 807.87(e)	DS-2.15.2.0-A-EN Safety and Setup Guide
Provide users with a surgical manual along with the instructions for use	N/A - The 3Shape Abutment Designer™ Software is not used in surgery
Provide all relevant precautions and warnings in the professional labelling	DS-2.15.2.0-A-EN Safety and Setup Guide
Precautions or warnings that relate to unpackaging or sterility	N/A - The 3Shape Abutment Designer™ Software is a software device and not supplied sterile.
If any parts are provided non-sterile we recommend that you provide sterilization instructions	N/A - The 3Shape Abutment Designer™ Software is a software device and not supplied sterile.
If patient labelling is appropriate, we recommend that you follow Guidance on Medical Device Patient Labelling; Final Guidance for Industry and FDA Reviewers	N/A - 3Shape Abutment Designer™ Software is a prescription device and patient labeling is not required.



3shape 

# Dental System™

Industry-leading scanning and CAD solutions



Excerpt  
version

# Get better esthetics and save costs with customized abutments

More and more labs are customizing abutments because they get improved clinical results, better esthetics and cost savings. Abutment Designer™ lets you automatically design the customized abutment and emergence profile with smooth transitions and optimal esthetics.




**All types of abutments**  
Create customized abutments, screw-retained crowns and anatomical abutments.

**COMING Auto-Abutment**  
Automatically generate an abutment that optimally fits the crown design.

**Easily adjust to Gingiva**  
Snap emergence profile to gingiva or use freeform tools.

**COMING Bone information**  
Improve abutment design by visualizing the jawbone.


**Abutment bevel**  
Easily add abutment bevel to your design for optimal crown support.



**NEW Add abutment position guides**  
Help your dentist easily and accurately place abutments, temporaries or veneers using custom designed positioning guides created in the same workflow.



**More than 50 implant libraries**  
Get abutments from both original and compatible-with implant manufacturers covering all global and key regional players.  
[www.3shapedental.com/implant-systems](http://www.3shapedental.com/implant-systems)



**COMING Bone information for better clinical results**  
Include bone data from CBCT scans to create better abutments without collisions and optimally shaped to the thickness of the gingiva.

## When accuracy and esthetics matter

Expand your lab's business by offering a high value service that is rapidly growing in demand. Design sophisticated implant bars and bridges with the utmost precision using flexible tools and ISO-documented 3Shape scanner accuracy.

**Design any type of bar**  
Support for a wide range of standard bar and bridge types such as Primary, Dolder, Hader, Hybrid, Canada, and Wrap-around.

**Extensive range of implant libraries\***  
Support for direct-seating and multi-unit abutments.

**Proven bar accuracy with 3Shape scanners**  
Accuracy on all 3Shape scanners is verified and documented according to ISO12836 and specially designed implant objects.

**NEW Freely apply different types of bars**  
Combine multiple types in a single design with the new smooth surface transition that ensures an optimal finish.

**NEW Improved finish with new Cut-to-gingiva**  
Match the bar to the underlying gingiva and achieve high quality with the newly enhanced surface finish.

**Add any type of attachment**  
Freely add slide or ball attachments, locators and retention holes. Adjust position and angle if needed.


\*Only libraries using FDA-cleared components permitted in the US.

**COMING Perfectly matching secondary structures**  
Design your bar and then easily and precisely create its matching secondary according to the bar's attachments and form.

**Esthetic implant bridges with gingiva**  
Easily design implant bridges complete with gingiva based on the final anatomy and a gingiva boundary that you define.

**COMING Designs obtained directly from the Denture or Wax-up**  
Scan the denture and virtually cut back to facilitate design of an optimal implant bridge, or scan an original wax-up bridge to create a file for copy milling.

## 3Shape develops 3D technologies for dental practices and labs



---

## Backing labs with care, technology and expertise



Expert  
Training



Unlimited  
Upgrade



Global  
Support

3Shape LABcare™ is an integral part of your annual 3Shape subscription. It bundles services that are designed to ensure your investment, secure maximum uptime, and help you get the most from your solution year after year.

Find 3Shape online



Questions? Contact FDA/CDRH/OCE/DID at [CDRH-FOISTATUS@fda.hhs.gov](mailto:CDRH-FOISTATUS@fda.hhs.gov) or 301.796.8118

DEPARTMENT OF HEALTH AND HUMAN SERVICES  
Food and Drug Administration

## Indications for Use

Form Approved: OMB No. 0910-0120  
Expiration Date: January 31, 2017  
See PRA Statement below.

510(k) Number (if known)  
K151455

Device Name  
3Shape Abutment Designer™ Software

### Indications for Use (Describe)

The 3Shape Abutment Designer Software is intended as an aid the restoration of chewing function in partially or fully edentulous mandibles and maxillae. The 3Shape Abutment Designer Software is intended for use by a dental practitioner or dental laboratory staff for designing the patient specific component of a two-piece, one-piece, or hybrid dental implant abutment. The single or multi-unit abutment design is intended to be used by the manufacturer of an endosseous dental implant abutment to create the final device.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## Software Requirement Specification

### 1. Introduction

Please refer to the Software Requirement Specifications enclosed in this volume.

Please note that the Validation Protocol templates are included in the Requirement Specifications.

For links between the Requirements Specifications and the Software Specifications, please see "VOL\_008\_Traceability Analysis" enclosed in this submission.

**Note:** Some of these documents may have requirements referring to other areas than 3Shape Abutment Designer™. See the Trace matrix for reference to which requirements are relevant for 3Shape Abutment Designer™.

### 2. Table of Contents

(b)(4)	Order Management (added document)
	Configuration
	Sending and receiving
	Manufacturing Interfaces
	Order Import and Export (added document)
	User Manual
	Material import export
	User Interface (newly created document)
	Parametric Abutment Design (updated version)
	Parametric Abutment Design. Version from 2.15.3.0 patch
	Wax-up Design
	Implant Bars (updated version)

Requirement Specification - (b)(4)

(b)(4)



(b)(4)





(b)(4)



(b)(4)





Requirement Specification - (b)(4)

(b)(4)

A large black rectangular redaction box covers the entire central portion of the page, obscuring all text and graphics that would otherwise be present.

(b)(4)



(b)(4)



Requirement Specification – (b)(4)

(b)(4)



(b)(4)



(b)(4)





(b)(4)



(b)(4)



(b)(4)



Requirement Specification - (b)(4)

(b)(4)



(b)(4)



(b)(4)




(b)(4)



(b)(4)







Requirement Specification - (b)(4)

(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)





Requirement Specification - (b)(4)

(b)(4)

A large black rectangular redaction box covers the entire central portion of the page, obscuring all text and graphics that would otherwise be present.



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)





3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

## Verification and Validation Documentation

(b)(4)





3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

(b)(4)





3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

(b)(4)





Validation Result – (b)(4) 

(b)(4)



(b)(4)




(b)(4)



(b)(4)





(b)(4)



Validation Result – (b)(4) 

(b)(4)



(b)(4)





(b)(4)



Requirement Specification – (b)(4)

(b)(4)



(b)(4)



(b)(4)





Validation Result – (b)(4)

(b)(4)

A large black rectangular redaction box covers the entire central portion of the page, obscuring all text and graphics that would otherwise be present.

(b)(4)



(b)(4)



(b)(4)





(b)(4)




(b)(4)



(b)(4)





(b)(4)

(b)(4)



(b)(4)



(b)(4)



(b)(4)





(b)(4)



(b)(4)



(b)(4)



(b)(4)





Validation Result – (b)(4)

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the section header and extending nearly to the bottom of the page.

(b)(4)



(b)(4)



(b)(4)





(b)(4)



(b)(4)



(b)(4)



(b)(4)




(b)(4)



Validation Result – (b)(4)

(b)(4)



(b)(4)



(b)(4)





(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)




(b)(4)



Verification Result – (b)(4)

(b)(4)



(b)(4)





(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)





(b)(4)



(b)(4)



(b)(4)





Verification Result – (b)(4)

(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)







Verification Result – (b)(4)

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the "Verification Result" header and extending nearly to the bottom of the page.

(b)(4)



(b)(4)



Verification Result – (b)(4)

(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)





(b)(4)



(b)(4)



(b)(4)





Verification Result – (b)(4)

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the section header and extending nearly to the bottom of the page. The text "(b)(4)" is written in red at the top left corner of this redacted area.

(b)(4)



(b)(4)



(b)(4)



(b)(4)





(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



Verification Result – (b)(4)

(b)(4)





(b)(4)



(b)(4)



(b)(4)



(b)(4)




(b)(4)



(b)(4)




Verification Result – (b)(4) 

(b)(4) 

(b)(4)







(b)(4)



---


Verification Protocol- (b)(4)

(b)(4)



(b)(4)





(b)(4)

(b)(4)



(b)(4)



(b)(4)





(b)(4)




(b)(4)





Verification Result - (b)(4)

(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)





(b)(4)



(b)(4)



Veirfication Result – (b)(4)

(b)(4)



(b)(4)



(b)(4)





Software Specification – (b)(4)

(b)(4)



(b)(4)



(b)(4)





(b)(4)



(b)(4)



(b)(4)



**PRODUCT RISK ASSESSMENT**

(b)(4)



(b)(4)



(b)(4)



(b)(4)





3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

---

## Handling of FDA Cleared Implant Libraries

### 1. Introduction

Please refer to the Standard Operating Procedure enclosed in this volume.

### 2. Table of Contents

(b)(4)	Handling of FDA Cleared Implant Libraries in 3Shape Software
--------	--





Standard Operating Procedure – (b)(4)

(b)(4)

A large black rectangular redaction box covers the majority of the page, obscuring all text and graphics that would otherwise be present in the main body of the document.

(b)(4)



















































































































3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION

---

## Appendix 2 – User Manual

### 1. Introduction

Please refer to the labelling enclosed in this appendix.

### 2. Table of Contents

Dental System Technical Guidelines
------------------------------------

(b)(4) - User Manual - Abutment Designer Except version
---



# Dental System Technical Guidelines

(b)(4) Draft Technical Guidelines



















# 3Shape Dental System 2015

## User Manual

(b)(4) Draft User Manual









































































































































## Appendix G: Implant Libraries for 3Shape Dental System™

Providers of Implant Solutions	Implant Systems						Implant Solutions				Manufacturing		Consumables (for local use)				Additional info	
	Nobel Biocare	Biomet3i	Straumann	Zimmer Dental	Astra Tech DENTSPLY	Friadent DENTSPLY CAMLOG	Other systems	1-Piece (Titanium)	1-Piece (Zirconia)	2-Piece (Ti base+ ZrAbut™)	Implant Bars & Bridges	Digital implant models	Providers milling center	Local milling by 3rd Party	Blanks with pre-milled interface	Titanium bases		Intraoral Scan bodies
<b>Alpha Bio-Tec. (IL)</b> www.alpha-bio.net							- Alpha Bio Tec	X	X	X	X		X	X	X			
<b>Argen (US)</b> www.argen.com	+	+	+	+	+	+		X	X			X		X	X	X	X	- FDA, CE
<b>ATLANTIS™ (SE)</b> www.dentsplyimplants.com	+	+	+	+	+	+	- Keystone Dental - BioHorizons	X	X			X						- FDA, CE
<b>BEGO (DE)</b> www.bego.com	+	+	+	+	+	+	- BEGO Implant Systems - and more	X	X	X	X	X			X	X		- FDA, CE - CE-labelled prosthetic screws
<b>BioComp (NL)</b> www.biocomp.eu							- BioComp	X	X				X	X	X	X	X	- CE, ISO
<b>Biodenta (CH)</b> www.biodenta.com	+	+	+	+	+	+	- Biodenta	X	X	X	X	X				X	X	- FDA, CE pending
<b>BioHorizons (US)</b> www.biohorizons.com	+		+	+			- BioHorizon	X	X				X	X	X	X		- FDA, CE
<b>Biomet (US)</b> www.biomet.com	+	+	+					X	X	X		X				X		- FDA, CE - Encode® healing abutments
<b>Biotech International (FR)</b> www.biotech-dental.com							- Biotech		X		X		X		X	X	X	- CE



<b>Digital Dental Group-DDG (ISR)</b> www.ddg-scanlab.com	+	+	+	+	+	+	+	+	+	- Microdent - BioHorizons - MIS-Implants - and more	X	X	X	X			X	X	X	X	X	- FDA pending, CE	
<b>DIO-Implants (KR)</b> www.dioimplant.com	+									- DIO Implants	X	X	X	X		X	X					- FDA, CE, KFDA, SFDA	
<b>Elos Medtech Pinol (DK)</b> www.elosmedtech.com	+	+	+	+	+	+	+	+	+	- Neoss	X		X	X	X	X	X			X	X	- FDA, CE	
<b>EuroTeknika (FR)</b> www.euroteknika.com	+	+	+	+	+					- Anthogyr - Biotech - Easy Implant - MIS Implants - and more	X		X	X	X	X			X	X	X	X	- FDA, CE
<b>GC Advanced Technologies (US)</b> www.gc-at.com	+	+	+	+	+					- BioHorizons - Sybron	X	X	X			X							
<b>Glidewell (US)</b> www.glidewell.com	+	+	+	+	+	+	+	+	+	- Prisma - DentalCraft - Keystone Dental - Neoss	X	X	X	X	X	X				X	X	X	- FDA, CE pending
<b>Heraeus Kulzer (DE)</b> www.heraeus.com	+	+	+	+	+	+	+	+	+	- Thommen	X	X	X		X	X							
<b>Ivoclar Wieland (DE)</b>	+	+	+	+	+	+	+	+	+	- Dentaurum Implants	X		X	X		X	X			X	X		
<b>LaStruttura (IT)</b> www.lastruttura.it	+	+	+	+	+	+	+	+	+	- Sweden&Martina - Megagen - Prodent - and many more	X	X		X	X	X					X	X	- CE 93/42
<b>Medentika (DE)</b> www.medentika.de	+	+	+	+	+	+	+	+	+	- Medentika M-Implant	X		X	X		X	X	X	X	X	X		- CE
<b>Medentis Medical (DE)</b> www.medentis.de	+		+	+	+					- ICX Templant	X		X	X	X	X				X	X		- CE
<b>Medical Production (FR)</b> www.medical-production.eu								+		- Euroteknika	X		X	X	X		X	X	X	X			- CEO 499 - ISO 9001 - ISO 13485
<b>MIS Implants Technology (IL)</b> www.mis-implants.com										- MIS Implants			X			X		X					- FDA, CE

<b>Neodent (BR)</b> www.neodent.com.br	+	+						- Neodent Implant Sy stems	X	X	X	X					X		X		
<b>Neoss (UK)</b> www.neoss.com	+	+	+	+	+			- Neoss	X		X	X	X	X	X		X	X	X	X	- FDA, CE
<b>Nobel Biocare (CH)</b> www.nobelbiocare.com	+	+	+	+	+				X	X	X			X							- FDA, CE
<b>NT Trading (DE)</b> www.nt-trading.com	+	+	+	+	+	+	+	- Thommen - Sweden&M artina - BEGO	X		X	X		X	X	X	X	X	X	X	- FDA, CE, CMDCAS, GOST R, TGC - Mexico Certification
<b>Phibo (ES)</b> www.phibo.com	+	+	+	+	+	+	+	- PHIBO - BTI - Sweden&M artina - MIS Implants - and many more	X	X	X	X	X	X					X	X	- FDA, CE - ISO 13485 - ISO 9001
<b>Prowital (DE)</b> www.prowital.de								- Prowital			X			X		X	X	X			- CE - EN-ISO 9001- V472008 - EN ISO- 13485- 2003AC2009 - Certificate Directive-93- 42 EWG
<b>Ritter Implants (DE)</b> www.ritterimplants.com			+	+				- Alpha- BioTec. - MIS Implants	X		X		X		X		X	X	X		- FDA, CE
<b>Straumann (CH)</b> www.straumann.com			+						X	X	X			X	X		X	X			- CE
<b>Sweden&amp;Martina S.p.A (IT)</b> www.swedenmartina.com								- Sweden&M artina	X	X	X	X	X	X			X	X	X		
<b>Target3D (USA)</b> www.target3d.com	+	+	+	+	+	+	+	- Osstem - Bio Horizons - MIS- Implants - and more	X	X	X	X	X		X	X	X	X	X		- FDA, CE
<b>Thommen Medical (CH)</b> www.thommenmedical.com								- Thommen			X	X					X	X			- FDA, CE
<b>TRI Implants (CH)</b> www.tri-implants.com								- TRI Implants			X						X	X	X		
<b>3dental (ES)</b> www.3dental.es	+	+	+	+	+	+		- BioComp - Sweden&M artina - DYNA	X	X		X	X	X						X	

+ -Available for Dental System™

x -Available from the manufacturer



**Note!** All information is given without guarantee and based exclusively on information made available by the implant system providers. Please contact your local 3Shape re-seller to obtain the most current list.



**Note!** Only Implant Libraries based on Dental Implants with 510(k) clearance can be used in Abutment Designer™ in the United States. Cleared libraries must be activated by 3Shape. Please ask your Implant Library Provider to contact 3Shape if you wish to use a cleared library that has not already been activated. The software will block any attempt to use a non-cleared library.



**Note!** In the United States, implant cases may only be manufactured by manufacturers that hold an implant abutment 510(k) or are milling per the specific instructions provided by the holder of a 510(k) for a patient specific implant abutment.

## Appendix H: Terms and Abbreviations

The following table explains file type abbreviations used in 3Shape Dental System.

File Format	Description
3ML	Zipped, compressed XML file typically used for setup and customization. 3Shape proprietary format.
3SE	3Shape and Sirona proprietary format for exported orders. Orders from Sirona system, can be imported into 3Shape system.
3OX	3Shape Communicate Order Exchange format used by partners and customers to retrieve and send orders. 3Shape proprietary format.
3OXZ	Zipped archive containing 3OX file and references to DCM models. 3Shape proprietary format.
3SI	3Shape and Sirona proprietary format for imported orders. Orders exported from 3Shape system, can be imported into Sirona system.
DCM	Dental Compressed Model file. Contains compressed 3D model data, attached objects (splines, annotations, etc.), marks and additional string properties. 3Shape proprietary format used for scans and CAD designs.
DLL	Dynamic-link library, Microsoft shared library concept.
DME	Dental System Material Export file. Contains materials, references to materials and external files. The file can be imported into another 3Shape Dental System. 3Shape proprietary format.
STL	Describes surface geometry of three-dimensional objects. Used for scans and CAD designs. Industry standard.
ULDC	3M Lava proprietary format. Order files from 3M scanners can be imported into 3Shape System.
XML	Extensible Markup Language, used for configuration files, etc.

## Appendix I : Contact Information

<b>3Shape Headquarters</b>	<b>3Shape North America</b>	<b>3Shape (Shanghai) Co., Ltd</b>	<b>3Shape South America</b>
Europe, Middle East & Africa Sales Holmens Kanal 7 1060 Copenhagen K Denmark	North American Sales Somerset Hills Corporate Center 10 Independence Boulevard, Suite 150 Warren, New Jersey 07059, USA	Asian Sales Room 906, Tower A of Eton Place No. 69, Dongfang Road 200120 Shanghai, China	Latin American and Caribbean Sales Carrera 13 # 82-91 Oficina 401 110221 Bogotá, Colombia
<b>P: +45 70 27 26 20</b>	<b>P: +1 908 867 0144</b>	<b>P: +86 21 5835 2281</b>	<b>P: +57 1 691 95 08</b>



3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

## **510(K) SUMMARY – Traditional 510(K)**

### **Submitter Information**

A Company Name: 3Shape A/S

B Company Address: Holmens Kanal 7  
DK-1060 Copenhagen K

C Company Phone: +45 7027 2620  
Company Fax: +45 7027 2621

D Contact Person: Hanne Nielsen  
Regulatory Affairs Manager

E Date Summary Prepared: June 28, 2016

### **Device Identification**

A Trade/proprietary Name: 3Shape Abutment Designer™  
Software

B Common Name: Abutment design software for dental  
laboratory

C Device Classification Name: Endosseous Dental Implant Abutment

C Regulation Number: 872.3630

C Classification: Class II

D Product Code: PNP

### **Predicate Device**

Sirona Dental CAD/CAM System (K100152).

### **Intended Use**

The 3Shape Abutment Designer Software is intended as an aid to the restoration of chewing function in partially or fully edentulous mandibles and maxillae. The 3Shape Abutment Designer Software is intended for use by a dental practitioner or dental laboratory staff for designing the patient specific component of a two-piece, one-piece or hybrid dental

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

implant abutment. The single or multi-unit abutment design is intended to be used by the manufacturer of an endosseous dental implant abutment to create the final device.

### **Device Description**

The 3Shape Abutment Designer™ Software receives scan data containing topographical characteristics of real teeth, dental impressions, or stone models. The software provides the user with the ability to create matching endosseous dental implant abutments using computer aided design.

The output of the device is a computer file containing the abutment(s) in digital form which can be used by manufacturers that hold an implant abutment 510(k) or are milling per the specific instructions provided by the holder of a 510(k) marketing clearance for a patient specific implant abutment. The 3Shape Abutment Designer Software output is restricted in the U.S. to be manufactured by a holder of an implant abutment 510(k) or milled per the specific instructions provided by the holder of a 510(k) marketing clearance for a patient specific implant abutment.

The 3Shape Abutment Designer™ Software can be run on properly configured "off-the-shelf" PC hardware running Microsoft® Windows and using a standard keyboard and mouse.

The 3Shape Abutment Designer™ Software requires the loading of implant libraries, which includes information such as implant type, maximum and minimum dimensional parameters for abutments, etc., created by separate abutment manufacturers and cleared by the FDA. In the US, the Implant Libraries are obtained via a 3Shape server after demonstration to 3Shape of the FDA clearance of the Implant Library.

### **Summary of the technological characteristics**

The 3Shape Abutment Designer™ is a software only device programmed in the Delphi programming language and has the following PC/laptop requirements and other technological characteristics as compared to the predicate device:



## 3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

Item	Submission Device Minimum Requirements	Submission Device Recommended	Predicate Device
OS	Windows 7 32-bit Professional*	Windows 7 64-bit Professional	Windows Vista 32-bit
RAM	4GB	8GB (16GB)	6GB
Video Card	512MB DirectX 10 (1GB DirectX 10)	1GB DirectX 11 (2GB DirectX 11)	512 MB DirectX 10 NVIDIA Quadro
Available HDD Space	250GB	500GB (1TB)	500GB
CPU	Intel Core i5 or equivalent	Intel Core i7 or equivalent	Intel Core i7
Monitor resolution	1440 x 900 pixels	1920 x 1080 pixels	Unknown
3D Mouse	None	3DConnexion SpaceMouse™	Unknown
Network	Internet connection		Unknown
Ports	USB 2.0 port for 3Shape desktop scanner		Serial Port PCI card for inEOS scanner
Mouse	Mouse with wheel button support		Unknown

## 3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

<b>Indications for Use</b>	The 3Shape Abutment Designer Software is intended as an aid in the restoration of chewing function in partially or fully edentulous mandibles and maxillae. The 3Shape Abutment Designer Software is intended for use by a dental practitioner or dental laboratory staff for designing the patient specific component of a two-piece, one-piece or hybrid dental implant abutment. The single or multi-unit abutment design is intended to be used by the manufacturer of an endosseous dental implant abutment to create the final device.	The Sirona Dental CAD/CAM System is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multi-unit cement retained restorations. The system consists of three major parts: TiBase, InCoris mesostructure, and CAD/CAM software. Specifically, the InCoris mesostructure and TiBase components make up a two-piece abutment which is used in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity. The InCoris mesostructure may also be used in conjunction with the CAMlog Titanium base CAD/CAM (types K2244.xxxx) (K083496) in the Camlog Implant System. The CAD/CAM software is intended to design and fabricate the InCoris mesostructure. The InCoris mesostructure and TiBase two-piece abutment is compatible with the following implant systems: Nobel Biocare Replace (K020646), Nobel Biocare Branemark (K022562), Friadent Xive (K013867), Biomet 3i Osseotite (K980549), Astra Tech Osseospeed (K091239), Zimmer Tapered Screw-Vent (K061410), and Straumann SynOcta (K061176).
<b>Software Output</b>	Digital encrypted or non-encrypted proprietary or .STL file including only the patient-specific abutment component for one-piece, two-piece, or hybrid abutment designs. The digital output does not include the abutment-to-implant connection interface.	.STL file of the ceramic mesostructure sent to Sirona Dental CAD/CAM System milling unit
<b>Physical - Output</b>	N/A – Submission device relies on separate regulatory clearance and manufacture of the abutment by a separate company	Two-piece TiBase abutment – pre-milled titanium base combined with ceramic mesostructure designed in Sirona CAD/CAM software.
<b>Milling Location</b>	Abutment Manufacturer or Dental laboratory per the 510(k) clearance of the dental abutment	Local milling of the ceramic abutment component.

The predicate Sirona device includes in the marketing clearance the two-piece TiBase Abutments (titanium bases and ceramic blocks) as a physical output as well as a validated milling unit and directions for assembly of the final dental abutment. The 3Shape Abutment Designer™ Software does not provide any physical parts that can come into contact with the patient. The 3Shape Abutment Designer Software instead provides only the digital design as an accessory to the physical dental abutment systems cleared by other manufacturers. The differences between the Indications for Use

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

Statement of the submission device and the predicate are related to the lack of physical output from the submission device and reliance on previous or subsequent FDA clearance of the physical abutment by a separate manufacturer. However, the difference does not change the intended use of the device.

### **Nonclinical Testing**

Software, hardware, and integration verification and validation testing was performed in accordance with the FDA Guidance Document "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices" (Issued on May 11, 2005) as well as the FDA Guidance Document "Off-The-Shelf Software Use in Medical Devices (Issued on September 9, 1999).

Prior to release, verification and validation testing of the 3Shape Abutment Designer™ Software has been completed using the approved acceptance criteria: Each user need has its own validation acceptance criteria; each specification has its own verification acceptance criteria; bug verification consists in ensuring issue is not reproducible; issues reported by beta partners must be reviewed and handled appropriately; beta partners must report functionality to be an improvement over previous version.

The validation suite includes validation of implemented mitigations related to device hazards identified in the risk management procedures.

### **Conclusion**

Based on a comparison of intended use, principle of operations, features and technical data, and the verification/validation test results, the 3Shape Abutment Designer™ Software is found to be substantially equivalent with the Predicate Device.

# Accuracy of 3Shape lab scanners

## Contents

Introduction.....	3
What is scanner accuracy? .....	3
Accuracy in dental applications.....	4
Determination of accuracy .....	4
Reference objects .....	4
Definition .....	4
Calibration and traceability .....	4
3Shape Calibration Objects .....	5
Gauge blocks.....	5
Reference sphere object.....	5
Implant bar validation object .....	5
ISO12836 Annex A – Inlay shaped specimen.....	6
ISO12836 Annex B – Crown and bridge shaped specimen.....	6
Free-form reference objects.....	6
Software .....	6
Production software.....	7
End user software.....	7
3 <sup>rd</sup> party software .....	7
3Shape lab scanner accuracy.....	8
Test program .....	8
Results .....	8
Implant bar case .....	8
ISO12836;2012, Annex A – Inlay Shape Specimen .....	11
ISO12836;2012, Annex B – Crown And Bridge Specimen .....	11
Appendix: Best practice for accurate scanning with 3Shape Lab Scanners .....	13
Introduction.....	13
Prerequisites.....	13

Accuracy of 3Shape lab scanners  
Issued: 15-09-2015



Preparations .....	13
Scanner off during nightttime.....	13
Scanner on during nightttime .....	13
Scanner calibration.....	13

## Introduction

### What is scanner accuracy?

The term *accuracy* is used in a number of senses, depending on the context and application. In the context of 3D scanning, accuracy relates to the capability of a 3D scanner to reproduce digitally the physical surface of a free form object.

In order to quantify this key property, the following measures are defined:

**Reproducibility:** The ability of a 3D scanner to reproduce a specific measurable dimension of a reference object of known dimensions. Mathematically, reproducibility is expressed as the standard deviation of a series of measurements of the dimension in question.

**Repeatability:** The ability of a 3D scanner to reproduce a specific measurable dimension of a reference object. Mathematically, repeatability is expressed as the difference between the average of a series of measurements (scans) and a true measured value. Repeatability is sometimes referred to as *trueness*.

See also ISO12836:2012: Dentistry – Digitizing devices for CAD/CAM systems for indirect dental restorations – Test methods for assessing accuracy

As for any physical quantity, dimensional measurements are affected by the environmental conditions in which they take place. For 3D scanning applications, key environmental conditions include relative humidity and temperature. All materials contract and expand with temperature; some more than others. The *coefficient of thermal expansion* is an intrinsic material property that dictates the sensitivity of a material to temperature changes.

However, not only will the object under investigation change size with temperature – also the 3D scanner system is subject to thermal effects; not least due to internal heat sources. In order to demonstrate the effect, consider the following example:

A caliper, made of steel, is used to measure the distance between two molars on a dental model in a laboratory. The caliper was produced at 20 degrees, and it is assumed that the scale on the ruler was absolutely correct at the time of production. In the laboratory, the temperature is now 25 degrees, and the dental model is produced at that temperature. The model is assumed to be an absolutely correct model of the patient teeth.

Using the caliper to measure the molar distance, the technician obtains the result 55 mm. However, steel has an expansion coefficient of around 12  $\mu\text{m}/\text{m}/\text{degree}$ , so in reality, the molar to molar distance has to be corrected by the amount of

$$5 \text{ deg} \times 12 \mu\text{m}/\text{m}/\text{deg} \times 0,055 \text{ m} = 3.3 \mu\text{m}$$

As temperature increases, the caliper expands. As the scale expands with the caliper, this means that any object measured with the caliper seems to be smaller than it really is; in case of decreasing temperature, an object will seem to be larger than it really is.

If the object was realized at a different temperature, this would either reduce or increase the magnitude of the induced measurement error, due to the dimensional change of the object itself.

So, the origin of the error from thermal effects is a combination of the object and the scanning system.

### **Accuracy in dental applications**

While the above example is trivial in some sense, it clearly demonstrates that effects due to temperature changes can easily be on the same order of magnitude as the accuracy specified by many 3D scanner system vendors. Therefore, any claim of scanner system accuracy must be with reference to test conditions in which the claimed accuracy can be obtained; including the allowed temperature interval from e.g. a system calibration temperature. With no such specifications, accuracy claims have no relevance to reality and should generally not be trusted.

As such effects scale linearly with the size of the case at hand, they may safely be ignored for small cases. However, in cases where large distances are involved, e.g. multi-unit bridges or implant bars, the end result will depend on the accuracy of the scanning and design process.

## **Determination of accuracy**

### **Reference objects**

#### **Definition**

As a reference object for 3D scanning may serve any object that represents the relevant scanning application. As such, in order to be relevant, the object must expose certain geometrical features that relate to the application. This can be free-form as well as standard geometrical shapes such as spheres, cylinders, cones etc.

#### **Calibration and traceability**

Most importantly, it shall be possible to determine the position, size and orientation of any such scannable feature using an application-independent measurement device. Frequently, coordinate measurement machines (CMMs) are used to perform such measurements. However, only to the extent that this reference measurement is traceable to national standards should it be used to provide objective evidence for the accuracy of the scanner in question. Hence, for any reference object, a calibration certificate referring to the accreditation of the measurement laboratory, the equipment used, and its calibration status should be provided.

(b)(4)





(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)





## Appendix 1 – Library Confirmation Letters

### 1. Introduction

Please refer to the Library Confirmation Letters enclosed in this appendix.

### 2. Table of Contents

Name	Comment
3ImplantSignature1	
3ShapeImplantSystem	Internal 3Shape library
ArgenImplantSystems	
AtlantisSystems	
BEGOSystemsFDA	
Bicon Implant Systems	
BiodentaCorporationImplants(C)FDA	
BiodentaCorporationImplants(E)	
BioHorizons	
BioHorizonsLibrary	
BlueSkyBioImplantSystems	
Camlog Implant Systems	
Camlog	
CAP Implant Systems	510k exempt
CAP Implants	510k exempt
CMC Implants	
Core3d Abutments Implant Library	
Core3d Implant Library Open	

Name	Comment
Core3d Implant Library	
Cortex Dental Implant System	
Creodental Implant System	New Library. Implant Library Clearance Document attached as well.
Dentium Implant System	New Library. Implant Library Clearance Document attached as well.
ETK Implant Library	
GlidewellAbutmentKit	Error in confirmation letter. K08340 should be K083480 and 141923 should be K141923
GlidewellAnatomyLibrary	Error in confirmation letter. K08340 should be K083480 and 141923 should be K141923
GlidewellBarKit	Error in confirmation letter. K08340 should be K083480 and 141923 should be K141923
GlidewellHybridKit	Error in confirmation letter. K08340 should be K083480 and 141923 should be K141923
Hiossen Implant Library	
ImplantDirectSystems	
Intra-lock implant library	
Medentika Implant Systems	
MegaGEN Implant System	New Library. Implant Library Clearance Document attached as well.
MIS Implant System	
Neodent Implant Systems	New Library. Implant Library Clearance Document attached as well.
Neoss Implant Library	
NobelBiocare Abutments	
NT-Trading Implants	
Paltop Implant Systems (C)	New Library. Implant Library Clearance Document attached as well.
Straumann One Piece Abutments	

Name	Comment
Sweden-Martina implant library	
TFI Systems Implant Library	New Library. Implant Library Clearance Document attached as well.
ThommenMedical Abutment kit	
TruCrown Implant System	New Library. Implant Library Clearance Document attached as well.
Zimmer Zfx Implant System	
Zimmer Zfx	
ZimmerAbutmentKit	
ZimmerImplantSignature1	
ZimmerImplantSystemC	
Zirlux Implant Systems	Error in confirmation letter. K08340 should be K083480 and 141923 should be K141923. New Library. Implant Library Clearance Document attached as well.

---

**Architecture Design Chart**

(b)(4)



TITLE: ARCHITECTURE DESIGN CHART

(b)(4)



TITLE: ARCHITECTURE DESIGN CHART

(b)(4)



TITLE: ARCHITECTURE DESIGN CHART

(b)(4)



TITLE: ARCHITECTURE DESIGN CHART

(b)(4)





TITLE: ARCHITECTURE DESIGN CHART

(b)(4)



May 9, 2016

## Creating an implant library for the 3Shape Dental System™



Version 1.0



May 9, 2016

Contents

Creating an implant library for the 3Shape Dental System™ ..... 1

1 Introduction ..... 4

2 3Shape’s implant library distribution policy ..... 4

3 Implant library and implant system concept ..... 4

4 The implant library installation workflow ..... 5

    4.1 Manual import of the implant library ..... 6

    4.2 Download over ftp ..... 7

5 Preparing CAD files for implant systems ..... 9

    5.1 Implant system parts ..... 9

        5.1.1 One-piece abutments ..... 9

        5.1.2 Two-piece abutments ..... 13

    5.2 3Shape Global Coordinate System for implant systems ..... 16

    5.3 Alignment of implant kit parts ..... 20

    5.4 Adjust number of triangles ..... 23

    5.5 Preparation of the base file ..... 25

    5.6 Preparation of the scan abutment ..... 31

6 Defining an Implant System in the Dental System™ ..... 34

    6.1 Importing CAD files and creating abutment kits ..... 35

    6.2 Rotate and translate a group of CAD files ..... 37

    6.3 Align scan abutment on top surface ..... 38

    6.4 Configuration of manufacturing options ..... 39

    6.5 Design matrix restrictions for abutments ..... 41

    6.6 Define scan abutment alignment points ..... 44

    6.7 Settings for implant analog friction bars ..... 46

    6.8 Implant system categories ..... 46

7 3Shape Global Implant Connection ID ..... 47

8 Implant library encryption ..... 50

**Available to all 3Shape users** ..... 50

**Available to selected 3Shape reseller** ..... 51

**Open CAD Output** ..... 51



May 9, 2016

<b>Encrypted CAD output</b> .....	51
9 Exporting an implant library.....	52
10 Creating Implant Libraries for use the United States.....	54
11 Dental System™ CAM Output .....	55
12 How to design a good scan abutment .....	56
13 How to see if a system is digitally signed .....	57
14 How to see if a library is FDA cleared.....	59
15 Example of Implant library .....	60

May 9, 2016

## 1 Introduction

In order to create an implant library for the design of customized abutments, implant bars and bridges in the Dental System™, 3Shape specific requirements must be met. This document is a detailed description of the requirements and will guide the implant library provider through the library creation.

## 2 3Shape's implant library distribution policy

3Shape is supporting the creation, distribution and promotion of original and compatible-with implant libraries under the following requirements:

- The implant library provider is the legal entity providing the Implant library.
- The implant library provider is responsible for obtaining all regulatory clearances relevant to the library and the parts thereof. Documentation for said clearance must be supplied to 3Shape where applicable.
- The implant library provider is responsible for the creation, the validation and any IP violation in relation to the library
- 3Shape only supports the distribution of the library to end customers if the providers has signed the implant library contract with 3Shape and the library has been provided to 3Shape
- The 3Shape reseller must approve the distribution of the implant library to his end-customers

## 3 Implant library and implant system concept

Implant Systems in the Dental System Control Panel are organized as following:

- An **implant library** is a digital representation of one or more **implant systems**.
- An **implant system** typically corresponds to a given implant type. In addition a number of implant systems can be grouped in a category, see sect. 5.5.
- Each **implant system** contains a number of **parts** – implants, scan abutments, screws, bases, interfaces, analog interfaces, etc. See sect. 4.1 for details.
- An **abutment kit** typically corresponds to a specific implant connection/diameter of the implant type. An abutment kit is an assembly of parts needed to complete the design of the customized abutment. Abutment kits are created by referencing the models added to a given implant system. Multiple abutment kits can reference the same parts and these kits, hence, allow you to easily combine the various parts belonging to a given implant system.

May 9, 2016

The diagram below illustrates the example of an implant system and one of its abutment kits:





Figure 3-1: Implant system structure in the Dental System™.

#### 4 The implant library installation workflow

Below figure illustrates the workflow for the installation of the implant libraries by the dental laboratory and the creation and provision of the implant library by the implant library provider.

May 9, 2016



Please see section 6 for a description how to create and update the implant library.

The dental lab has two ways to retrieve the implant library from the implant library provider:

1. Manual import of the implant library
2. Download over ftp

#### 4.1 Manual import of the implant library

In the Control Panel the under Tools → Import / export it is possible to import the implant library by selecting “Import materials”.




Figure 4-1: Implant library import in the Dental System™.

May 9, 2016

After locating and selecting the implant library file a dialogue shows the contents i.e. the implant systems that the implant library file contains. The appropriate implant systems are then selected for import.




Figure 4-2: Import and update materials.

## 4.2 Download over ftp

In the Control Panel the under Tools → Download center it is possible to retrieve materials from ftp either by selecting the Download materials or the Download libraries option.



May 9, 2016




Figure 4-3: Materials download.




Figure 4-4: Implant library download.

In the Download libraries option it is possible to select which library should be retrieved and installed. Please see section 14 on how to identify if a library is cleared for use in the US by FDA.



May 9, 2016

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the date and extending nearly to the bottom of the page.

May 9, 2016



### **Implant**

There are no specific requirements for the model of the implant. This model is used for visualization purposes only.



### **Scan abutment**

The model corresponds to the physical scan abutment. However, the non-visible part of the model -such as the inside of the screw hole and the interface- should be deleted to maximize alignment accuracy.

The design of the scan abutment is important for achieving an accurate fit, which is particular important for implant bridges.




### **Screw**

There are no specific requirements for the model of the screw. The model of the screw is used to automatically generate a correctly shaped hole through the abutment.

In case of using a pre-milled blank with pre-defined screw hole, the generated screw hole can be automatically removed from the CAM output (see sect. 3.2).

(b)(4)



**Blank**

The blank model should be a cylinder without screw hole and interface. This model is used to visualize the size of the blank.




**Digital model analog interface (optional)**

Digital model analog interface is used for creating a hole, where model analog is inserted. The model should contain either

- **two** boundaries (model analog interface with the **bottom** insertion)

OR

- **one** boundary at the top of the model (model



May 9, 2016

analog interface with insertion from the **top**)



**Inner limit (optional)**

The inner limit model should be a cylinder without screw hole. This model is used to visualize the limitations for the minimum geometry of the abutment.

May 9, 2016

## 5.1.2 Two-piece abutments



### Implant

There are no specific requirements for the model of the implant. This model is used for visualization purposes only.



### Scan abutment

The model corresponds to the physical scan abutment. However, the non-visible part of the model -such as the inside of the screw hole and the interface - should be deleted to maximize alignment accuracy.

The design of the scan abutment is important for achieving an accurate fit, which is particular important for implant bridges.



### Screw

There are no specific requirements for the model of the screw. The model of the screw is used to automatically generate a correctly shaped hole through the abutment.

In case of using a pre-milled blank with pre-defined screw hole, the generated screw hole can be automatically removed from the CAM output (see sect. 3.2).

(b)(4)




#### **Abutment interface**

There are no specific requirements for the model of the interface.

*This file is for visualization only, the 3Shape software does not allow any changes by the end-user to this file during the abutment design process.*

May 9, 2016



**Blank (optional)**

The blank model should be a cylinder without screw hole and interface. This model is used to visualize the size of the blank.

**Digital model analog interface (optional)**


Digital model analog interface is used for creating a hole, model analog is inserted.

The model should contain either

- **two** boundaries (model analog interface with the **bottom** insertion)

OR

- **one** boundary at the top of the model (model analog interface with insertion from the **top**)



**Inner limit (optional)**

The inner limit model should be a cylinder without a screw hole. This model is used to visualize the limitations for the minimum geometry of the abutment.





(b)(4)





(b)(4)





(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the 3shape logo and extending nearly to the bottom of the page.



(b)(4)





(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the (b)(4) label and extending nearly to the bottom of the page.



(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the 3shape logo and extending nearly to the bottom of the page.

(b)(4)





(b)(4)







(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the 3shape logo and ending above the footer. The text "(b)(4)" is written in red at the top left corner of this redacted area.



(b)(4)





(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the (b)(4) label and extending nearly to the bottom of the page.



May 9, 2016

(b)(4)





(b)(4)





(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the 3shape logo and extending nearly to the bottom of the page.



(b)(4)





May 9, 2016

(b)(4)







May 9, 2016

(b)(4)





(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the 3shape logo and extending nearly to the bottom of the page.



May 9, 2016

## **6 Defining an Implant System in the Dental System™**

Creating implant systems is easily done using the Dental System™ Control Panel. How to configure abutment kits for the customized abutment design is described step by step in the following section.

Please note, that implant libraries are forward compatible.

If you create a library, using Dental System 2014 version, it won't be possible to use it in older version of the software – Dental System 2013, while Dental System 2014 and newer would be able to import it.

We are constantly adding new features, such as 3Shape Global Coordinate System for implant systems (see chapter 4.2), and 3Shape Global Implant Connection ID (see chapter 6), so it's important to use the latest version in order to use new functions. However, older version would provide implant library with wider compatibility over different versions.

May 9, 2016

## 6.1 Importing CAD files and creating abutment kits

Open the Dental System™ Control Panel and go to **Abutments > Implant Systems**. This will bring up the **Details** mode of the implant systems:




Figure 5-5: Defining implant systems and abutment kits in the Control Panel

1. To add a new implant system, press the **Add** button under Implant System - **(1)** - and enter a unique ID, the name for the new system and a version number.
2. To add parts/CAD files to the selected implant system, press **Add** under the implant system parts - **(3)** - and browse the file. Import all the desired parts before proceeding to the next step.

May 9, 2016

- To create abutment kits from the parts of the active implant system, click the **Add** button under Abutment kits - (2). This will bring up a **wizard** which will guide you through the part selection.

A **Table mode** for the implant systems has been implemented for a better overview and an easier and faster editing of abutment kits. A double-click activates the fields and values can be changed.




Figure 5-6: Table mode for implant systems


Browsing CAD files and creating abutment kits have still to be done in the details mode as shown and described above.

May 9, 2016

## 6.2 Rotate and translate a group of CAD files

New feature "Transform", allows rotating and repositioning of entire set of models with the same transformation saving time and reducing the number of errors.

Select folder with models, and click "Apply". Files would be updated in the same folder.



May 9, 2016

### 6.3 Align scan abutment on top surface

The alignment algorithm has been improved in order to achieve better results mainly for the alignment of scan abutments with a horizontal top surface. This improvement ensures a better vertical locking of the implant position as well as the locking of the scan abutment rotation.

Implant system providers can enable the option *Align with top surface* when browsing their scan abutment CAD file (see sect. 5.1 step 2).




Figure 5-7: Align with top surface option

To read more about the requirements of scan abutments in order to achieve optimized scan results, please refer to sect. 10.

May 9, 2016

## 6.4 Configuration of manufacturing options

The abutment manufacturing options allow the user to specify how abutments designed using the abutment kit in question should be exported.

The screenshot shows a dialog box titled "Edit abutment kit" with the following sections and controls:

- Identification:** Three text input fields for "ID", "Name", and "External ID".
- 3Shape Global Coordinate system:** A checkbox labeled "Support 3Shape Global Coordinate system" with a note: "Check this option if you want 3Shape Global Coordinate system to be supported."
- Abutment export options:**
  - Checked checkbox: "Use implant coordinates" with a note: "Check this option if you want the abutment to be saved in the implant coordinate system." Below it is a spinner control for "Rotate output around Z-axis" set to 0 degrees.
  - Unchecked checkbox: "Keep base" with a note: "Check this option if you want to keep the shape of the abutment base. Uncheck this option if you want to remove the base geometry from the exported abutment."
  - Checked checkbox: "Append hole patches" with a note: "Check this option if you want to remove the screw hole from the exported abutment."
- Advanced options:** Three spinner controls for "Vertical screw offset", "Extra drill hole radius", and "Hole fillet", all set to 0 mm. A note for "Hole fillet" reads: "Avoid thin and sharp walls around screw hole exit".
- CAD block options:** Two spinner controls for "Z axis offset" and "Rotate around implant", both set to 0 mm and 0 deg respectively.

At the bottom are "Cancel", "Back", and "Next" buttons.

Figure 5-10: Manufacturing options

**Use implant coordinates** - when checked, the final abutment geometry is stored in the common coordinate system of the implant system parts. Furthermore, it allows rotating the output around the z-axis, if required for the production. If this option is not checked, the final implant will be stored in the coordinate system of the preparation scan used for designing the abutment.




May 9, 2016


**Keep base** - when unchecked, the base geometry is removed from the exported model and replaced with a planar surface.

**Note:** For regulatory reasons, the 3Shape software will refuse to output the connector base if the end user is located in the United States.

This allows the use of milling blanks with the pre-milled base geometry. The base geometry will be removed for one-piece abutments and kept for two-piece abutments.




Base geometry removed




Base geometry included

**Append hole patches** - when checked, the screw hole of the abutment is closed in the exported geometry. This allows the use of milling blanks with the pre-drilled screw holes.



Screw hole closed



Screw hole appended



May 9, 2016

(b)(4)

A large black rectangular redaction box covers the majority of the page content. The text "(b)(4)" is written in red at the top left corner of this box.



May 9, 2016

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the date and ending above the footer. The text "(b)(4)" is written in red at the top left corner of this redacted area.



May 9, 2016

(b)(4)



May 9, 2016

## 6.6 Define scan abutment alignment points

Implant system providers can now pre-define alignment points to ensure an optimal alignment of the scan abutment during the scanning process.

The alignment points can be set under Preview (1) in the Abutment Kit box. In the preview window, 1-point or 3-point alignment (2) can be selected. By clicking on *Set point* or *Set points*, the alignment points can be pre-defined on the scan abutment.





Figure 5-13: Set alignment points


The following should be taken into account when pre-defining the scan abutment alignment points:

- The alignment points can be set on any sufficient and visible surface – the placement in deep indentations as e.g. screw channels are to be avoided
- The alignment points should not be placed close to or on any edge
- The placement of the points should be considered as an easy-to-match position for the user

May 9, 2016



1-point alignment




3-point alignment


Figure 5-14: Recommended alignment point setting

May 9, 2016

## 6.7 Settings for implant analog friction bars

When importing a model analog for the production of implant models additional parameters can be defined as shown below.



Details	
ID	15045_AnalogInterfacePW <span>Change ID</span>
Name	AnalogCAD
External ID	Not available
Digital signatures	<none>
CAD file	 C:\3Shape Configuration\DentalSystem\Dental System Control Panel\library\3Shape\Abutments\ImplantSystems\Impor
Friction surface top	0.0 <input type="text"/> mm
Friction surface bottom	0.0 <input type="text"/> mm


**Friction surface top** - The distance from Y= 0 mm to the top of the friction bar. If for instance you want the friction bar to have its top point in Y= -1 mm, set the parameter to 1.0 mm.

**Friction surface bottom** - The distance from Y= 0 mm to the bottom of the friction bar. If for instance you want the friction bar to have its bottom point in Y= -6 mm, set the parameter to 6.0 mm.

## 6.8 Implant system categories

Implant systems can be organized into categories for an easier selection in the Dental Manager Order form.

1. Click the **Add** button in the *Implant system categories* window.
2. Enter a name for your category in the appeared form.
3. Click **OK**.
4. Click the **Add current system** button to include the implant system you are currently working with into the selected category.
5. Every implant system page contains **Implant system categories**, so you can click on your category and add that implant system to the list.



May 9, 2016

## 7 3Shape Global Implant Connection ID

In order to ensure the automatic mapping of implant selections between the 3Shape applications Dental System™, TRIOS® and Implant Studio™, a Global Implant Connection ID has to be defined for each original implant brand.

Implant system providers are responsible for generating Global Implant Connection IDs when creating the implant library in the Dental System™ Control Panel.

By providing the requested information as described below, the Global Implant Connection ID will be automatically generated according to the naming convention of the 3Shape Global ID:

*Manufacturer\_System\_Connection*

Manufacturer	3Shape
System	Demo
Connection	4.3



**Global Implant Connection ID**  
**3Shape\_Demo\_4.3**



May 9, 2016




Figure 6-15: Creation of Global Implant Connection ID in the Control Panel

The following information has to be provided by the implant system provider:

- (1) Original Manufacturer:** Select one of the listed implant manufacturers (e.g. *Nobel Biocare, Straumann, Phibo or Avinent*)
- (2) System:** Insert the implant brand (e.g. *NobelActive, SLA, TSA Advance or Coral HI*)
- (3) Connection:** Type the name or dimension of the interface or platform on the implant shoulder level (e.g. *NP, RN, S4 or 4.1*)



May 9, 2016

Example of Global Implant Connection IDs:

Manufacturer	System	Connection	Global Implant Connection ID
Straumann	Soft Tissue Level	RN	Straumann_Soft Tissue Level_RN
Straumann	Soft Tissue Level	WN	Straumann_Soft Tissue Level_WN
Straumann	Soft Tissue Level	NNC	Straumann_Soft Tissue Level_NNC
Straumann	Bone Level	RC	Straumann_Bone Level_RC
Straumann	Bone Level	NC	Straumann_Bone Level_NC
Nobel Biocare	Replace	NP	Nobel Biocare_Replace_NP
Nobel Biocare	Replace	RP	Nobel Biocare_Replace_RP
Nobel Biocare	Replace	WP	Nobel Biocare_Replace_WP
Nobel Biocare	Replace	6.0	Nobel Biocare_Replace_6.0
Nobel Biocare	Branemark	NP	Nobel Biocare_Branemark_NP
Nobel Biocare	Branemark	RP	Nobel Biocare_Branemark_RP
Nobel Biocare	Branemark	WP	Nobel Biocare_Branemark_WP
Nobel Biocare	Conical Connection	NP	Nobel Biocare_Conical Connection_NP
Nobel Biocare	Conical Connection	RP	Nobel Biocare_Conical Connection_RP
Nobel Biocare	Conical Connection	30	Nobel Biocare_Conical Connection_3.0
BEGO	S-Line	5.5	BEGO_S-Line_5.5
Biodenta	BL	B0	Biodenta_BL_B0
CAMLOG	iSy	3.8	CAMLOG_iSy_3.8
Phibo	TSA/ADV	S5	Phibo_TSA/ADV_S5
Euroteknika	Naturactis	NP	Euroteknika_Naturactis_NP

**Important:** In order to ensure the compatibility of implant libraries with the 3Shape applications Dental System™, TRIOS® and Implant Studio™, the implant libraries must be created in conformity with the requirements of the 3Shape Global Coordinate System described in sect. 4.2.



May 9, 2016

(b)(4)





May 9, 2016

(b)(4)

A large, solid black rectangular redaction box covers the majority of the page's content. The text "(b)(4)" is written in red at the top left corner of this redacted area.

May 9, 2016

## 9 Exporting an implant library

The export functionality in Dental System™ Control Panel provides an easy exchange of materials. Go to **Tools > Export/Import** button and click on “Export materials”



Select the desired materials from the list. Before exporting the implant library the following options can be selected:

- Lock exported materials
- Always overwrite materials with the same IDs when importing
- Delete pre-existing items in destination materials when importing

**Lock exported materials-** when this option is checked, the material settings cannot be changed at the receiving end (recommended).

**FDA cleared libraries are locked:** All implant libraries that have been cleared for use in the US are automatically locked so that the end user is not able to change the implant library settings.

**Always overwrite materials with the same IDs when importing-** when this option is checked, materials with the same ID will automatically be replaced at the receiving end.

**Delete pre-existing items in destination materials-** when this option is checked, all existing materials at the receiving end will be deleted when importing new materials.

May 9, 2016

Click on the “Export to File” button to export the implant library in .dme format.





May 9, 2016

## **10 Creating Implant Libraries for use the United States**

If the end-user is located in the United States, only FDA-cleared abutment components/libraries are available for use in 3Shape Dental System™.

This is due to a request from the FDA and thus only applies for installations in the United States.

In order to make a library available in Dental System™ the Implant Library Provider is required to send documentation to 3Shape in form of an electronic copy of the relevant 510(k) clearance letters.

The Implant Library Provider also needs to sign and forward a confirmation letter stating that the library only consists of FDA-cleared components to the extent that clearance is required for said component.

Please be aware that the Implant Library Provider is responsible for ensuring that the library only contains FDA-cleared components to the extent that clearance is required for said component.

Implant Libraries for use in the United States will only function in the 3Shape Dental System™ if they are encrypted and locked, and have been approved by 3Shape.

Please contact [implant@3shape.com](mailto:implant@3shape.com) for more information.



May 9, 2016

(b)(4)





May 9, 2016

## 12 How to design a good scan abutment

When designing a scan body it is very important to understand how the scanners work. Basically the scan body is scanned and the scanned surface is matched to the CAD of the scan flag minimizing the difference between the 2 surfaces. During this matching it is critical that the scanned surface has sufficient surfaces to lock all positions and orientations.



Dental System™ 2012 - Model Builder™ (min. 2.14)

<http://www.youtube.com/watch?v=kQ8DA6wD9aM&list=UUxEI9LrUI3A7SjkOPbDLPpg&index=3&feature=plcp>

The following should take in account when designing a scan abutment for the 3Shape Dental System™:

- Sufficient horizontal surface for vertical locking of the implant position
- Conical scan abutment are to be avoid in order to achieve better results in the vertical alignment
- A non-symmetric feature, preferred a limited vertical flat surface or sufficient bevel surface to lock the rotation of the scan abutment
- Avoid sharp edges to ensure an accurate alignment of the scan abutment and the CAD file of the scan abutment
- Scanable material, e.g. Peek, which has sufficient scan properties, but notice that there are different types of peek with varying scan quality
- Make sure to use intraoral approved material to ensure the usability for intraoral and lab scanner

May 9, 2016

## 13 How to see if a system is digitally signed

Digital signature is essential for protection of intellectual rights and for FDA clearance mechanism.

Please ensure, that your implant library models are digitally signed with following steps:

1. Open Dental System Control Panel, and navigate to Abutments -> Implant Systems. Select your implant system.
2. Choose model file – implant model, for example, and scroll down to details.
3. Check “Digital signature” value.

Dental System Control Panel - Implant system

Home Back Save Help Table mode

Abutments > Implant systems

Implant system 1

3Shape Implant 4.3 ...

Edit Implant system

ID 3S\_ImplantSystem Change ID

Name 3Shape Implant 4.3 mm

Version 1.0

Abutment kits

2 pcs Abutment 4.3 mm  
1 pcs Abutment 4.3 mm  
Implant Bar/Bridge 4.3mm

Implant models 2

3Shape 4.3mm implant

Scan abutment models

3Shape scan abutment

Screw models

3Shape screw

scroll down

May 9, 2016

**Details**

ID	3S_ImplantSystem_Implant	<input type="button" value="Change ID"/>
Name	<input type="text" value="3Shape 4.3mm implant"/>	
External ID	<input type="text" value="Not available"/>	
Global Implant Connection ID	<input type="text"/>	
Original Manufacturer	<input type="text" value="Other..."/>	
System	<input type="text"/>	
Connection	<input type="text"/>	
Scan flag design	<input type="text"/>	
Color	<input type="text"/>	
Optional	<input type="text"/>	
Digital signatures	<input type="text" value="&lt;none&gt;"/> 3	

This is an example of a file that is **NOT signed**

**Details**

ID	3S_ImplantSystem_Implant	<input type="button" value="Change ID"/>
Name	<input type="text" value="3Shape 4.3mm implant"/>	
External ID	<input type="text" value="Not available"/>	
Global Implant Connection ID	<input type="text"/>	
Original Manufacturer	<input type="text" value="Other..."/>	
System	<input type="text"/>	
Connection	<input type="text"/>	
Scan flag design	<input type="text"/>	
Color	<input type="text"/>	
Optional	<input type="text"/>	
Digital signatures	<input type="text" value="3ShapeImplantSystem"/> 3	

**Digitally signed.** The signature used is "3ShapeImplantSystem"

May 9, 2016

## 14 How to see if a library is FDA cleared

1. Open Dental System Control Panel, and navigate to Abutments -> Implant Systems. Select your implant system.
2. Choose abutment kit and scroll down to details.
3. Check FDA logo presence.

Dental System Control Panel - Implant system

Home Back Save Help Table mode

Abutments > Implant systems

Implant system 1

3Shape Implant 4.3 ...

Recycle bin

Details

Total: 211 items

Edit Implant system

ID: 3S\_ImplantSystem Change ID

Name: 3Shape Implant 4.3 mm

Version: 1.0

Abutment kits 2

2 pcs Abutment 4.3 mm

1 pcs Abutment 4.3 mm

Implant Bar/Bridge 4.3mm

scroll down

Implant models

3Shape 4.3mm implant


Scan abutment models

3Shape scan abutment

Screw models

3Shape screw

May 9, 2016

<b>Details</b> FDA 3	
ID	3S_AbutmentKit_2pc_4_3 <input type="button" value="Change ID"/>
Name	2 pcs Abutment 4.3 mm
External ID	50000
Implant ID	3S_ImplantSystem_Implant 

Please note that the library in the figures above is the 3Shape internal test and demo library which cannot be used for production.

## 15 Example of Implant library

This section provides an example of an Implant library setting; the standard 3Shape Implant 4.3mm library which is provided with the software.

The library contains two abutment kits:

- 1 pcs Abutment 4.3 mm
- 2 pcs Abutment 4.3 mm

By selecting on of the kits the contents of the implant kit is shown. Below illustrates the 1 pcs Abutment 4.3 mm implant kit content and preview of the 3d models:

May 9, 2016

Dental System Control Panel - Implant system

Home Back Save Help Table mode Transform

Abutments > Implant systems

**Implant system**

3Shape Implant 4.3 mm

+ Add  
 + Copy  
 X Delete  
 ↑ Move up  
 ↓ Move down

**Recycle bin**

No items  
 Restore all  
 Empty the Recycle Bin

**Details**

Total: 5 items

**Edit Implant system**

ID: 3S\_ImplantSystem Change ID

Name: 3Shape Implant 4.3 mm

Version: 1.0

**Abutment kits**

1 pcs Abutment 4.3 mm  
Implant Bar/Bridge 4.3mm  
2 pcs Abutment 4.3 mm

+ Add  
 + Copy  
 X Delete  
 Edit  
 Preview

**Implant models**

3Shape 4.3mm implant

+ Add  
X Delete

**Scan abutment models**

3Shape scan abutment

+ Add  
X Delete

**Screw models**

3Shape screw  
30124002\_screw

+ Add  
X Delete

**Abutment base models**

3Shape 4.3 mm 2 pc...  
3Shape 4.3 mm 1 pc...  
3Shape 4.3 mm bar

+ Add  
X Delete

**Abutment interface models**

3Shape 4.3mm interface

+ Add  
X Delete

**Analog interface models**

3Shape analog interface

+ Add  
X Delete

**Blanks**

3Shape Blank 13x12 mm

+ Add  
X Delete

**Inner limits**

+ Add  
X Delete

**Anchor pins**

+ Add  
X Delete

**Details** *FDA*

ID: 3S\_AbutmentKit\_1pc\_4\_3 Change ID

Name: 1 pcs Abutment 4.3 mm

External ID:

Implant ID: 3S\_ImplantSystem\_Implant

Scan abutment ID: 3S\_ImplantSystem\_ScanFlag

Scan abutment display name: 3Shape scan abutment

Align with top surface:

Screw ID: 3S\_ImplantSystem\_Screw


Abutment base ID: 3S\_ImplantSystem\_Base1pcs

Analog interface ID: 3S\_ImplantSystem\_Analog

Blank ID: 3S\_ImplantSystem\_Blank

Anchor pin ID:

May 9, 2016



Below shows the settings for the 3Shape 1 pcs Abutment 4.3 mm implant kit. For details on the specific settings and how to maintain these please see section 6.

May 9, 2016

Edit abutment kit

×

### Identification

ID	<input type="text" value="3S_AbutmentKit_1pc_4_3"/>
Name	<input type="text" value="1 pcs Abutment 4.3 mm"/>
External ID	<input type="text"/>

### 3Shape Global Coordinate system

Support 3Shape Global Coordinate system  
Check this option if you want 3Shape Global Coordinate system to be supported.

### Abutment export options

Use implant coordinates  
Check this option if you want the abutment to be saved in the implant coordinate system.

Rotate output around Z-axis  deg

Keep base  
Check this option if you want to keep the shape of the abutment base  
Uncheck this option if you want to remove the base geometry from the exported abutment.

Append hole patches  
Check this option if you want to remove the screw hole from the exported abutment.

### Advanced options

Vertical screw offset	<input type="text" value="0"/> mm
Extra drill hole radius	<input type="text" value="0"/> mm
Hole fillet <small>Avoid thin and sharp walls around screw hole exit</small>	<input type="text" value="0"/> mm

### CAD block options

Z-axis offset	<input type="text" value="0"/> mm
Rotate around implant	<input type="text" value="0"/> deg

Cancel

Back

Next



May 9, 2016

Edit abutment kit



**Implants**

3Shape 4.3mm implant

ID 3S\_ImplantSystem\_Implant

Name 3Shape 4.3mm implant

CAD model



Cancel

Back

Next

May 9, 2016

Edit abutment kit

×

**Scan Abutments**

3Shape scan abutment

Please select the different scan abutments which this kit can use

ID 3S\_ImplantSystem\_ScanFlag

Name 3Shape scan abutment

CAD model



Cancel

Back

Next

May 9, 2016

Edit abutment kit

×

**Screws**

3Shape screw  
30124002\_screw

ID 3S\_ImplantSystem\_Screw

Name 3Shape screw

CAD model



Cancel

Back

Next

May 9, 2016

Edit abutment kit



**Screws**

- 3Shape screw
- 30124002\_screw**

ID 43000\_Screw4277

Name 30124002\_screw

CAD model



Cancel

Back

Next

May 9, 2016

Edit abutment kit



**Interfaces (optional)**

3Shape 4.3mm interface

ID 3S\_ImplantSystem\_Interface

Name 3Shape 4.3mm interface

CAD model



Cancel

Back

Next


May 9, 2016

Edit abutment kit



**Analog interfaces (optional)**

3Shape analog interface



ID 3S\_ImplantSystem\_Analog

Name 3Shape analog interface

CAD model



Cancel

Back

Next

May 9, 2016

Edit abutment kit



**Inner limits (optional)**



Cancel

Back

Next

May 9, 2016

Edit abutment kit

×

**Blanks (optional)**

3Shape Blank 13x12 mm

ID 3S\_ImplantSystem\_Blank

Name 3Shape Blank 13x12 mm

CAD model



Cancel

Back

Next



May 9, 2016

Edit abutment kit



**Anchor pins (optional)**



Cancel

Back

Next

May 9, 2016

Edit abutment kit

×

**Materials**

- Zirkon
- Ti
- Wax
- PMMA
- CoCr
- Model Material
- DREVE FotoDentModel Beige
- 3D Systems Visijet
- Prefab crown
- priticrown
- Titan
- Zirkon (Configuration\_01)

Please select the different materials in which this kit can be manufactured

Cancel

Back

Next



May 9, 2016

Edit abutment kit

×

**Bases**

- 3Shape 4.3 mm 2 pcs abutm...
- 3Shape 4.3 mm 1 pcs abutm...
- 3Shape 4.3 mm bar

ID 3S\_ImplantSystem\_Base1pcs  
Name 3Shape 4.3 mm 1 pcs abutment  
CAD model



Cancel

Back

Next

# 3Shape Dental System 2015

## User Manual

(b)(4) Draft User Manual









































































































































## Appendix G Implan Libraries for 3Shape Denta System™

The latest list of Implan Libraries for 3Shape Denta System can be found on this web page:  
[http://support.3shape.com/media/547752/Implan\\_libraries\\_for\\_the\\_3shape\\_denta\\_system.pdf](http://support.3shape.com/media/547752/Implan_libraries_for_the_3shape_denta_system.pdf).



**Note!** All information is given without guarantee and based exclusively on information made available by the implant system provider. Please contact your local 3Shape reseller to obtain the most current list.



**Note!** Only Implant Libraries based on Dental Implants with 510(k) clearance can be used in Abutment Design r™ in the United States. Cleared libraries must be activated by 3Shape. Please ask your Implant Library Provider to contact 3Shape if you wish to use a cleared library that is not already been activated. The software will block any attempt to use a non-cleared library.



**Warning** In the United States, important cases may not be manufactured by manufacturers that do not have an important abutment 51(k) or are dental laboratories milling per the specific instructions for a patient specific important abutment provided by the holder of a 51(k) specifically cleared for dental laboratory milling. Please check with the implant library provider to ensure an appropriate location for sending the digital file.

## Appendix H: Terms and Abbreviations

The following table explains file type abbreviations used in 3Shape Dental System.

File Format	Description
3ML	Zipped, compressed XML file typically used for setup and customization. 3Shape proprietary format.
3SE	3Shape and Sirona proprietary format for exported orders. Orders from Sirona system, can be imported into 3Shape system.
3OX	3Shape Communicate Order Exchange format used by partners and customers to retrieve and send orders. 3Shape proprietary format.
3OXZ	Zipped archive containing 3OX file and references to DCM models. 3Shape proprietary format.
3SI	3Shape and Sirona proprietary format for imported orders. Orders exported from 3Shape system, can be imported into Sirona system.
DCM	Dental Compressed Model file. Contains compressed 3D model data, attached objects (splines, annotations, etc.), marks and additional string properties. 3Shape proprietary format used for scans and CAD designs.
DLL	Dynamic-link library, Microsoft shared library concept.
DME	Dental System Material Export file. Contains materials, references to materials and external files. The file can be imported into another 3Shape Dental System. 3Shape proprietary format.
STL	Describes surface geometry of three-dimensional objects. Used for scans and CAD designs. Industry standard.
ULDC	3M Lava proprietary format. Order files from 3M scanners can be imported into 3Shape System.
XML	Extensible Markup Language, used for configuration files, etc.

## Appendix I : Contact Information

<b>3Shape Headquarters</b>	<b>3Shape North America</b>	<b>3Shape (Shanghai) Co., Ltd</b>	<b>3Shape South America</b>
Europe, Middle East & Africa Sales Holmens Kanal 7 1060 Copenhagen K Denmark	North American Sales Somerset Hills Corporate Center 10 Independence Boulevard, Suite 150 Warren, New Jersey 07059, USA	Asian Sales Room 906, Tower A of Eton Place No. 69, Dongfang Road 200120 Shanghai, China	Latin American and Caribbean Sales Carrera 13 # 82-91 Oficina 401 110221 Bogotá, Colombia
<b>P: +45 70 27 26 20</b>	<b>P: +1 908 867 0144</b>	<b>P: +86 21 5835 2281</b>	<b>P: +57 1 691 95 08</b>



(b)(4)



(b)(4)





# Implant libraries for the 3Shape Dental System™

May, 2016

Implant systems	Implant systems							Implant Solutions					Manufacturing		Consumables (for local use)					Additional info				
	Nobel Biocare	Biomet 3i	Straumann	Zimmer Dental	DENTSPLY - Astra Tech	DENTSPLY - Friadent	CAMLOG	Other systems	1- Piece (Titanium)	1- Piece (Zirconia)	2- Piece (Ti base + Zr abutment)	Implant Bars & Bridges	Digital Implant models	Providers milling center	Local milling by 3rd party	Blanks with pre-milled interface (for 1-piece abutment)	Titanium bases (for 2-piece abutments)	Intraoral Scan bodies	Digital model analogs		Multi-unit abutments	Cleared for use in US	Comments	
3dental (ES) www.3dental.es	3shape reseller	👍	👍	👍	👍	👍	👍	- BioComp - Sweden & Martina - DYNA	X	X		X	X	X					X					
AB-Dental (IL) www.ab-dent.com				👍				- MIS Alpha Bio	X		X		X			X	X	X	X	X		- CE		
Abutments4life (DE) www.abutments4life.de	3shape reseller	👍	👍	👍	👍	👍	👍	- BioHorizons - Osstem - ICX-Templant - and more	X		X	X	X	X	X	X	X	X	X			- CE		
Alpha-Bio Tec (IL) www.alpha-bio.net								- Alpha Bio Tec	X		X	X	X		X	X						- CE		
Argen (US) www.argen.com	3shape reseller	👍	👍	👍	👍	👍	👍		X		X			X		X	X	X	X		X	- CE		
ATLANTIS™ (SE) - FDA cleared only www.dentsplyimplants.com	3shape reseller	Awaiting information from the provider					👍	👍	Awaiting information from the provider													X	- CE	
ATLANTIS™ (SE) - all libraries www.dentsplyimplants.com	3shape reseller	👍	👍	👍	👍	👍	👍	- DENTSPLY ANKYLOS XIVE® - FRALIT™ - Keystone Dental - BioHorizons	X	X				X								- CE		
BEGO (DE) www.bego.com	3shape reseller	👍	👍	👍	👍	👍	👍	- BEGO Implant Systems - and more	X	X	X	X		X			X	X				- CE - CE-labelled prosthetic screws		
BioComp (NL) www.biocomp.eu								- BioComp	X		X				X	X	X	X	X			- CE ISO		
Bicon (US) www.bicon.com								- Bicon			X	X			X			X			X	- CE Health Canada		
Biodenta (CH) - FDA cleared only www.biodenta.com	3shape reseller	Awaiting information from the provider							- Biodenta	Awaiting information from the provider													X	- CE pending
Biodenta (CH) - all libraries www.biodenta.com	3shape reseller	👍	👍	👍	👍	👍	👍	- Biodenta	X		X	X	X	X				X	X			- CE pending		
Biomet (US) - FDA cleared only www.biomet.com	3shape reseller	Awaiting information from the provider		👍				Awaiting information from the provider													X	- CE - encode healing abutments		
Biomet (US) - all libraries www.biomet.com	3shape reseller	👍	👍	👍					X	X	X			X				X				- CE - encode healing abutments		
Biotech Dental (FR) www.biotech-dental.com								- Biotech			X		X		X		X	X	X	X		- CE		
BlueSkyBio (US) blueskybio.com		👍		👍	👍	👍		- BlueSkyBio	X		X				X	X	X			X	X	- CE, Health Canada		
Bredent (DE) www.bredent.com								- Sky			X				X		X	X	X			- CE		
CADBLU (US) www.cadbludental.com	3shape reseller	👍	👍	👍	👍	👍		- Bio Horizons - MIS Implants	X		X	X	X	X				X	X					
CADstar (AT) www.cadstar.dental		👍	👍	👍	👍	👍	👍	- BREIDENT DENTSPLY NIVE - CONOLOG IMPLANTS ICK - NEOS PROWITAL SIC - ALPHATECH and others	X		X	X		X			X		X	X		- CE (0297) - ISO 13485 - Certificate Directive 93-42-EEC - (Annex II)		
CAMLOG (CH) - FDA cleared only www.camlog.com						👍		- CONLOG - Isy	Awaiting information from the provider													X		
CAMLOG (CH) - all libraries www.camlog.com						👍		- CONLOG - Isy	X		X	X		X	X		X	X						
CAP (US) www.cap-us.com	3shape reseller	👍	👍	👍	👍	👍			X		X	X	X	X	X	X	X	X	X		X			

**3Shape Headquarters**  
Europe, Middle East & Africa  
Holmens Kanal 7  
1060 Copenhagen, Denmark  
Tel: +45 7027 2620

**3Shape Asia**  
Room 906, Tower A of Eton Place  
No. 69, Dongfang Road  
200120 Shanghai, China  
Tel: +86 21 5835 2281

**3Shape Latin America**  
Carrera 13 # 82-91  
Oficina 401  
110221 Bogotá, Colombia  
Tel: +57 1691 9508

**3Shape North America**  
Somerset Hills Corporate Center  
10 Independence Boulevard, Suite 150  
Warren, New Jersey 07059, USA  
Tel: +1 (908) 867 0144

**More information at**  
www.3shape.com  
www.3shapedental.com  
info@3shape.com



Implant systems	Implant systems							Implant Solutions					Manu- facturing		Consumables (for local use)					Additional info			
	Nobel Biocare	Biomet 3i	Straumann	Zimmer Dental	DENTSPLY - Astra Tech	DENTSPLY - Friadent	CAMLOG	Other systems	1- Piece (Titanium)	1- Piece (Zirconia)	2- Piece (Ti base+ Zr abutment)	Implant Bars & Bridges	Digital Implant models	Providers milling center	Local milling by 3rd party	Blanks with pre-milled interface (for 1-piece abutment)	Titanium bases (for 2-piece abutments)	Intraoral Scan bodies	Digital model analogs		Multi-unit abutments	Cleared for use in US	Comments
<b>CMC (US)</b> - FDA cleared only www.custom-milling.com	👍	👍	👍	👍	👍	👍	👍	- BioHorizons - Osstem - LifeCore - Intra-Lock - Neoss	X	X	X	X		X				X			X	- CE	
<b>CMC (US)</b> - all libraries www.custom-milling.com	👍	👍	👍	👍	👍	👍	👍	- BioHorizons - Osstem - LifeCore - Intra-Lock - Neoss	X	X	X	X		X				X				- CE	
<b>Conmet (RU)</b> conmet.ru								- Conmet	X		X		X			X	X	X				- Dental implants have a bioactive surface - CE ISO 9001 and ISO 13485	
<b>Conexão Sistemas de Protese (BR)</b> www.conexao.com.br	👍		👍					- Conexão Sistema de Protese	X	X	X			X		X	X	X	X			- ANVISA CE	
<b>Core3D International (NL)</b> - FDA cleared only www.core3dcenters.com	Awaiting information from the provider							Awaiting information from the provider														X	- CE Health Canada TGA (AUS) Taiwan regulatory - MHLW:Ministry of Health Labour and Welfare (JP)
<b>Core3D International (NL)</b> - all libraries www.core3dcenters.com	👍	👍	👍	👍	👍	👍	👍	- Avinent - BioComp - BioHorizons - MIS Implants - and many more	X		X	X	X	X	X	X	X	X	X	X		- CE Health Canada TGA (AUS) Taiwan regulatory - MHLW:Ministry of Health Labour and Welfare (JP)	
<b>Cortex Dental (IL)</b> www.cortex-dental.com				👍				- Cortex Dental			X	X			X	X	X	X	X		X	- CE	
<b>Creodont (US)</b> www.Creomc.com	👍	👍	👍	👍	👍			- Hiossen	X		X	X		X		X					X		
<b>C-Tech Implant (IT)</b> www.c-tech-implant.com								- C-Tech Implant	X		X	X		X		X				X		- CE	
<b>Degudent (DE)</b> www.degudent.com	👍	👍	👍	👍	👍	👍	👍	- Medentika	X	X	X			X	X	X							
<b>Dental Consulting (DE)</b> www.gadau-consulting.com	👍		👍	👍	👍	👍	👍	- ICI-Medentis - Osstem	X		X	X	X	X	X	X	X	X	X			- CE	
<b>Dentaurum Implants (DE)</b> www.dentaurum-implants.de								- Dentaurum Implants			X				X	X	X					- CE PAL	
<b>Dentegris Deutschland (DE)</b> www.dentegris.de								- Dentegris Implants			X	X			X	X							
<b>Dentium (KR)</b> www.dentium.com								- Implantium SuperLine - SimpleLine II - NR Line Custom Abutment	X	X	X	X		X	X						X	- CE CFDA PMDA TFDA GOST	
<b>Dentsply-Friadent (DE)</b> www.dentsply-friadent.com						👍		- Xivo - Ankylos	X	X	X			X								- CE	
<b>DESS (ES)</b> www.dess-abutments.com	👍	👍	👍	👍	👍	👍		- BioHorizons - Ankylos - Mis	X		X	X		X	X	X	X			X		- CE - ISO 13485 - ISO 9001 - Health Canada pending.	
<b>DIO-Implants (KR)</b> www.dioimplant.com	👍							- DIO Implants	X	X	X	X		X	X							- CE, KFDA, SFDA	
<b>Dynamic Abutment Solutions (ES)</b> www.dynamicaabutment.com	👍	👍	👍	👍	👍	👍	👍	- Osstem - MIS Implants - Magagen - Biohorizons - and more	X		X	X		X	X	X	X			X		- CE, ISO 9001, ISO13485, CMDCAS	
<b>Easy Implant (FR)</b> www.easyimplant.com									X		X	X		X		X				X		- CE	
<b>Elos Medtech Pinol (DK)</b> www.elosmedtech.com	👍	👍	👍	👍	👍	👍	👍	- Neoss	X		X	X	X	X	X	X		X	X			- CE	
<b>Euroteknika (FR)</b> www.euroteknika.com	👍	👍	👍	👍	👍			- Anthogyr - Biotech - Easy Implant - MIS Implants - and more	X		X	X	X	X		X	X	X	X	X	X	- CE	
<b>GC Advanced Technologies (US)</b> www.gc-at.com	👍	👍	👍	👍	👍			- BioHorizons - Sybron	X	X	X			X									

**3Shape Headquarters**  
 Europe, Middle East & Africa  
 Holmens Kanal 7  
 1060 Copenhagen, Denmark  
 Tel: +45 7027 2620

**3Shape Asia**  
 Room 906, Tower A of Eton Place  
 No. 69, Dongfang Road  
 200120 Shanghai, China  
 Tel: +86 21 5835 2281

**3Shape Latin America**  
 Carrera 13 # 82-91  
 Oficina 401  
 110221 Bogotá, Colombia  
 Tel: +57 1691 9508

**3Shape North America**  
 Somerset Hills Corporate Center  
 10 Independence Boulevard, Suite 150  
 Warren, New Jersey 07059, USA  
 Tel: +1 (908) 867 0144

**More information at**  
[www.3shape.com](http://www.3shape.com)  
[www.3shapedental.com](http://www.3shapedental.com)  
[info@3shape.com](mailto:info@3shape.com)



Implant systems	Implant systems							Implant Solutions					Manu- facturing	Consumables (for local use)					Additional info					
	Nobel Biocare	Blomet 3i	Strueman	Zimmer Dental	DENTSPLY - Astra Tech	DENTSPLY - Friadent	CAMLOG	Other systems	1- Piece (Titanium)	1- Piece (Zirconia)	2- Piece (Ti base + Zr abutment)	Implant Bars & Bridges	Digital implant models	Providers milling center	Local milling by 3rd party	Blanks with pre-milled interface (for 1-piece abutment)	Titanium bases (for 2-piece abutments)	Intraoral Scan bodies	Digital model analogs	Multi-unit abutments	Cleared for use in US	Comments		
<b>Glidewell (US)</b> www.glidewelldental.com	👍	👍	👍	👍	👍	👍	👍	- Prisma 6 - DentalCraft - Keystone Dental - Neos	X	X	X	X	X	X			X	X	X		X	- CE pending		
<b>Green DenTech / Denrade (TW)</b> www.denrade.com	👍	👍	👍	👍	👍	👍		- Biohorizons - Dentium - Osstem and more...	X							X	X						- CE	
<b>Heraeus Kulzer (DE)</b> www.heraeus.com	👍	👍	👍	👍	👍	👍	👍	- Thommen	X	X	X		X	X										
<b>Implant Direct (US)</b> www.implantdirect.com/custom-direct	👍		👍	👍				- Biohorizons	X			X	X	X					X		X		- CE Health Canada	
<b>Intra-Lock International (US)</b> www.intra-lock.com								- INTRA-LOCK	X		X	X				X	X	X	X	X	X	X	- CE Health Canada	
<b>Ivoclar Wieland (DE)</b> www.wieland-dental.de	👍	👍	👍	👍	👍	👍		- Dentaurum Implants	X		X	X		X	X		X		X					
<b>Lasak (CZ)</b> www.lasak.com	👍		👍		👍			- Lasak BiomQ - Lasak Implant	X	X	X	X		X		X	X				X		- CE - ISO 13485 - ISO 9001	
<b>LaStruttura (IT)</b> www.lastruttura.it	👍	👍	👍	👍	👍	👍		- Sweden & Martina - Magagn - Prudent - and many more	X	X		X	X	X				X	X				- CE 93/42	
<b>Medentika (DE)</b> www.medentika.de	👍	👍	👍	👍	👍	👍	👍	- Medentika M- Implant	X		X	X		X	X	X	X	X			X		- CE	
<b>Medentis Medical (DE)</b> www.medentis.de	👍		👍		👍	👍		- ICX Templant	X		X	X	X	X			X	X					- CE	
<b>Medical Instinct (DE)</b> www.medical-instinct.de								- Medical Instinct	X	X	X	X			X	X	X	X			X		- CE	
<b>Médical Production (FR)</b> www.medical-production.eu		👍	👍	👍	👍	👍		- Eurotenika	X		X	X	X		X	X	X	X			X		- CE0499 - ISO 9001 - ISO 13485	
<b>Megagen Implant (KR)</b> www.imegagen.com								- AnyRidge Internal - AnyOne Internal - MINI Internal - Octa Level			X		X	X	X		X	X	X		X		- CE - KFDA	
<b>MIS Implants Technology (IL)</b> - FDA cleared only www.mis-implants.com								- MIS Implants	Awaiting information from the provider											X		- CE		
<b>MIS Implants Technology (IL)</b> - all libraries www.mis-implants.com								- MIS Implants	X		X			X	X	X	X	X		X			- CE	
<b>Neodent (BR)</b> www.neodent.com.br	👍		👍					- Neodent Implant Systems	X	X	X	X			X			X				X		
<b>Neoss (UK)</b> - FDA cleared only www.neoss.com	Awaiting information from the provider							- Neoss	Awaiting information from the provider											X		- CE		
<b>Neoss (UK)</b> - all libraries www.neoss.com	👍	👍	👍	👍	👍	👍		- Neoss	X		X	X	X	X	X	X	X	X	X	X	X	X	- CE	
<b>Nobel Biocare (CH)</b> www.nobelbiocare.com	👍	👍	👍	👍	👍				X	X	X			X								X	- CE	
<b>NT Trading (DE)</b> www.nt-trading.com	👍	👍	👍	👍	👍	👍	👍	- Thommen - Sweden & Martina - BEGO	X		X	X		X	X	X	X	X	X		X		- CE, CMCAS, GOST R, TGC - Mexico Certification	
<b>Phibo (ES)</b> www.phibo.com	👍	👍	👍	👍	👍	👍	👍	- BTI - Sweden & Martina - MIS Implants - and many more	X	X	X	X	X	X					X	X			- CE - ISO 13485 - ISO 9001	
<b>Prowital (DE)</b> www.prowital.de								- Prowital			X				X		X	X	X				- CE - EN-ISO 9001-V472008 - EN ISO-13485-2003AC2009 - Certificate Directive-93-42 IWG	
<b>Ritter Implants (DE)</b> www.ritterimplants.com			👍	👍				- Alpha -Bio Tec. - MIS Implants	X		X		X		X		X	X	X				- CE	

**3Shape Headquarters**  
Europe, Middle East & Africa  
Holmens Kanal 7  
1060 Copenhagen, Denmark  
Tel: +45 7027 2620

**3Shape Asia**  
Room 906, Tower A of Eton Place  
No. 69, Dongfang Road  
200120 Shanghai, China  
Tel: +86 21 5835 2281

**3Shape Latin America**  
Carrera 13 # 82-91  
Oficina 401  
110221 Bogotá, Colombia  
Tel: +57 1691 9508

**3Shape North America**  
Somerset Hills Corporate Center  
10 Independence Boulevard, Suite 150  
Warren, New Jersey 07059, USA  
Tel: +1 (908) 867 0144

**More information at**  
www.3shape.com  
www.3shapedental.com  
info@3shape.com





Implant systems	Implant systems							Implant Solutions					Manu- facturing		Consumables (for local use)					Additional info			
	Nobel Biocare	Biomet 3i	Straumann	Zimmer Dental	DENTSPLY - Astra Tech	DENTSPLY - Friadent	CAMLOG	Other systems	1- Piece (Titanium)	1- Piece (Zirconia)	2- Piece (Ti base + Zr abutment)	Implant Bars & Bridges	Digital Implant models	Providers milling center	Local milling by 3rd party	Blanks with pre-milled interface (for 1-piece abutment)	Titanium bases (for 2-piece abutments)	Intraoral Scan bodies	Digital model analogs		Multi-unit abutments	Cleared for use in US	Comments
<b>Providers of Implant solutions</b>																							
Straumann (CH) - FDA cleared only www.straumann.com			👍					Awaiting information from the provider														X	- CE
Straumann (CH) - all libraries www.straumann.com			👍					X	X	X				X	X		X	X	X				- CE
Sweden&Martina S.p.A (IT) www.sweden-martina.com							- Sweden&Martina	X	X	X	X	X	X				X	X	X				
Target3D (US) www.target3d.com	👍	👍	👍	👍	👍	👍	👍	- Osstem - BioHorizons - MIS-Implants - and more	X	X	X	X	X	X	X	X	X	X	X	X	X		- CE
T.F.I. System srl (IT) www.tfisystem.it							- Easy Grip			X	X				X		X	X	X	X	X	X	- CE
Thommen Medical (CH) www.thommenmedical.com							- Thommen			X	X						X	X				X	- CE
TRI Implants (CH) www.tri-implants.com							- TRI Implants			X							X	X	X				
Vulcan Custom Dental (US) www.vulcandental.com				👍			- BioHorizons	X		X	X	X	X	X	X	X	X	X	X	X	X	X	
Zfx (DE) www.zfx-dental.com	👍	👍	👍	👍	👍	👍	👍	- Megagen - Thommen Medical - Microdent - Osstem - and more	X	X	X	X	X	X	X		X			X			- CE
Zimmer Dental (US) www.zimmerdental.com				👍				X	X					X								X	- CE

**Please note:**

- The above information is solely provided by third parties, and 3Shape assumes no responsibility for its accuracy or completeness.

**3Shape Headquarters**  
Europe, Middle East & Africa  
Holmens Kanal 7  
1060 Copenhagen, Denmark  
Tel: +45 7027 2620

**3Shape Asia**  
Room 906, Tower A of Eton Place  
No. 69, Dongfang Road  
200120 Shanghai, China  
Tel: +86 21 5835 2281

**3Shape Latin America**  
Carrera 13 # 82-91  
Oficina 401  
110221 Bogotá, Colombia  
Tel: +57 1691 9508

**3Shape North America**  
Somerset Hills Corporate Center  
10 Independence Boulevard, Suite 150  
Warren, New Jersey 07059, USA  
Tel: +1 (908) 867 0144

**More information at**  
www.3shape.com  
www.3shapedental.com  
info@3shape.com



# Implant libraries for the 3Shape Dental System™

June, 2016

Implant systems	Implant systems							Implant Solutions					Manu- facturing		Consumables (for local use)					Additional info			
	Nobel Biocare	Biomet 3i	Straumann	Zimmer Dental	DENTSPLY - Astra Tech	DENTSPLY - Friadent	CAMLOG	Other systems	1- Piece (Titanium)	1- Piece (Zirconia)	2- Piece (Ti base + Zr abutment)	Implant Bars & Bridges	Digital Implant models	Providers milling center	Local milling by 3rd party	Blanks with pre-milled interface (for 1-piece abutment)	Titanium bases (for 2-piece abutments)	Intraoral Scan bodies	Digital model analogs		Multi-unit abutments	Cleared for use in US	Comments
3dental (ES) www.3dental.es	reseller	👍	👍	👍	👍	👍	👍	- BioComp - Sweden & Martina - DYNA	X	X		X	X	X					X				
AB-Dental (IL) www.ab-dent.com				👍				- MIS Alpha Bio	X		X		X			X	X	X	X	X		- CE	
Abutments4life (DE) www.abutments4life.de		👍	👍	👍	👍	👍	👍	- BioHorizons - Osstem - ICX-Templant - and more	X		X	X	X	X	X	X	X	X	X			- CE	
Alpha-Bio Tec (IL) www.alpha-bio.net								- Alpha Bio Tec.	X		X	X	X		X	X						- CE	
Argen (US) www.argen.com	reseller	👍	👍	👍	👍	👍	👍				X			X		X	X	X	X		X	- CE	
ATLANTIS™ (SE) - FDA cleared only www.dentsplyimplants.com	reseller	Awaiting information from the provider				👍	👍	Awaiting information from the provider														X	- CE
ATLANTIS™ (SE) - all libraries www.dentsplyimplants.com	reseller	👍	👍	👍	👍	👍	👍	- DENTSPLY ANKYLOS XIVE® - FRALIT™ - Keystone Dental - BioHorizons	X	X				X								- CE	
BEGO (DE) www.bego.com	reseller	👍	👍	👍	👍	👍	👍	- BEGO Implant Systems - and more	X	X	X	X		X			X	X				- CE - CE-labelled prosthetic screws	
BioComp (NL) www.biocomp.eu								- BioComp	X		X				X	X	X	X	X			- CE ISO	
Bicon (US) www.bicon.com								- Bicon			X	X			X			X			X	- CE Health Canada	
Biodenta (CH) - FDA cleared only www.biodenta.com	reseller	👍	👍	👍	👍	👍		- Straumann Bone Level (no Tissue Level) - BioHorizons Internal Osstem (Hiossen) TSH - Biodenta Bone Level and Tapered	X		X	X	X	X	X	X	X	X	X	X	X	- Implant Bridges only available for Abutment Level not Implant level. - 3rd party milling center must become Biodenta Contract Manufacturer and comply to QSR	
Biodenta (CH) - all libraries www.biodenta.com	reseller	👍	👍	👍	👍	👍		- Straumann Bone Level (no Tissue Level) - BioHorizons Internal Osstem (Hiossen) TSH - Biodenta Bone Level and Tapered	X		X	X	X	X	X	X	X	X	X	X	X	- CE CFDA TFDA	
Biomet (US) - FDA cleared only www.biomet.com	reseller	Awaiting information from the provider		👍				Awaiting information from the provider														X	- CE - encode healing abutments
Biomet (US) - all libraries www.biomet.com	reseller	👍	👍	👍					X	X	X			X				X				- CE - encode healing abutments	
Biotech Dental (FR) www.biotech-dental.com								- Biotech			X		X		X		X	X	X	X		- CE	
BlueSkyBio (US) blueskybio.com		👍		👍	👍	👍		- BlueSkyBio	X		X				X	X	X				X	X	- CE Health Canada
Bredent (DE) www.bredent.com								- Sky			X				X		X	X	X			- CE	
CADBLU (US) www.cadbludental.com	reseller	👍	👍	👍	👍	👍		- Bio Horizons - MIS Implants	X		X	X	X	X				X	X				
CADstar (AT) www.cadstar.dental		👍	👍	👍	👍	👍	👍	- BREDED DENTSPLY XIVE - CONELOG MEDENTS ICX - NESS PROWITAL SIC - ALPHATECH and others	X		X	X		X			X		X	X		- CE (0297) - ISO 13485 - Certificate Directive 93-42-EEC - (Annex II)	
CAMLOG (CH) - FDA cleared only www.camlog.com							👍	- CONELOG - Isy									X	X			X		
CAMLOG (CH) - all libraries www.camlog.com							👍	- CONELOG - Isy	X		X	X		X	X		X	X					

**3Shape Headquarters**  
Europe, Middle East & Africa  
Holmens Kanal 7  
1060 Copenhagen, Denmark  
Tel: +45 7027 2620

**3Shape Asia**  
Room 906, Tower A of Eton Place  
No. 69, Dongfang Road  
200120 Shanghai, China  
Tel: +86 21 5835 2281

**3Shape Latin America**  
Carrera 13 # 82-91  
Oficina 401  
110221 Bogotá, Colombia  
Tel: +57 1691 9508

**3Shape North America**  
Somerset Hills Corporate Center  
10 Independence Boulevard, Suite 150  
Warren, New Jersey 07059, USA  
Tel: +1 (908) 867 0144

**More informat on at**  
www.3shape.com  
www.3shapedental.com  
info@3shape.com



Implant systems	Implant systems							Implant Solutions					Manu- facturing		Consumables (for local use)					Additional info				
	Providers of implant solutions	Nobel Biocare	Biomet 3i	Straumann	Zimmer Dental	DENTSPLY - Astra Tech	DENTSPLY - Friadent	CAMLOG	Other systems	1- Piece (Titanium)	1- Piece (Zirconia)	2- Piece (Ti base + Zr abutment)	Implant Bars & Bridges	Digital implant models	Providers milling center	Local milling by 3rd party	Blanks with pre-milled interface (for 1-piece abutment)	Titanium bases (for 2-piece abutments)	Intraoral Scan bodies		Digital model analogs	Multi-unit abutments	Cleared for use in US	Comments
CAP (US) www.cap-us.com	reseller	👍	👍	👍	👍	👍	👍			X		X	X	X	X	X	X			X				
CMC (US) - FDA cleared only www.custom-milling.com	reseller	👍	👍	👍	👍	👍	👍		- BioHorizons - Osstem - LifeCore - Intra-Lock - Neoss	X	X	X	X		X				X			X	- CE	
CMC (US) - all libraries www.custom-milling.com	reseller	👍	👍	👍	👍	👍	👍		- BioHorizons - Osstem - LifeCore - Intra-Lock - Neoss	X	X	X	X		X				X				- CE	
Conmet (RU) conmet.ru									- Conmet	X		X		X					X	X	X		- Dental implants have a bioactive surface - CE ISO 9001 and ISO 13485	
Conexão Sistemas de Protese (BR) www.conexao.com.br		👍		👍					- Conexão Sistema de Protese	X	X	X			X				X	X	X	X	- ANVISA CE	
Core3D International (NL) - FDA cleared only www.core3dcenters.com	reseller	Awaiting information from the provider							Awaiting information from the provider														X	- CE Health Canada TGA (AUS) Taiwan regulatory - MHLW:Ministry of Health Labour and Welfare (JP)
Core3D International (NL) - all libraries www.core3dcenters.com	reseller	👍	👍	👍	👍	👍	👍	👍	- Avinent - BioComp - BioHorizons - MIS Implants - and many more	X		X	X	X	X	X	X	X	X	X	X		- CE Health Canada TGA (AUS) Taiwan regulatory - MHLW:Ministry of Health Labour and Welfare (JP)	
Cortex Dental (IL) www.cortex-dental.com					👍				- Cortex Dental			X	X			X	X	X	X	X		X	- CE	
Creodent (US) www.Creomc.com		👍	👍	👍	👍	👍			- Hiossen	X		X	X		X			X				X		
C-Tech Implant (IT) www.c-tech-implant.com									- C-Tech Implant	X		X	X		X			X				X	- CE	
Degudent (DE) www.degudent.com	reseller	👍	👍	👍	👍	👍	👍	👍	- Medentika	X	X	X			X	X		X						
Dental Consulting (DE) www.gadau-consulting.com		👍		👍	👍	👍	👍	👍	- ICK-Medents - Osstem	X		X	X	X	X	X	X	X	X	X			- CE	
Dentaurum Implants (DE) www.dentaurum-implants.de									- Dentaurum Implants			X				X		X	X				- CE PAL	
Dentegris Deutschland (DE) www.dentegris.de									- Dentegris Implants			X	X			X		X						
Dentium (KR) www.dentium.com									- Implantium - SuperLine - SimpleLine II - NR Line Custom Abutment	X	X	X	X		X	X						X	- CE CFDA PMDA TFDA GOST	
Dentsply-Friadent (DE) www.dentsplyfriadent.com	reseller					👍			- Xive - Ankylos	X	X	X			X								- CE	
DESS (ES) www.dess-abutments.com		👍	👍	👍	👍	👍	👍		- BioHorizons - Ankylos - MIs	X		X	X		X	X		X	X		X		- CE - ISO 13485 - ISO 9001 - Health Canada pending.	
DIO-Implants (KR) www.dioimplant.com	reseller	👍							- DIO Implants	X	X	X	X		X	X							- CE KFDA SFDA	
Dynamic Abutment Solutions (ES) www.dynamicaabutment.com		👍	👍	👍	👍	👍	👍	👍	- Osstem - MIS Implants - Megagen - Biohorizons - and more	X		X	X		X	X	X	X			X		- CE, ISO 9001, ISO13485,CMDCAS	
Easy Implant (FR) www.easyimplant.com										X		X			X			X				X	- CE	
Elos Medtech Pinol (DK) www.elosmedtech.com		👍	👍	👍	👍	👍	👍	👍	- Neoss	X		X	X	X	X	X	X		X	X			- CE	
Euroteknika (FR) www.euroteknika.com		👍	👍	👍	👍	👍			- Anthogyr - Biotech - Easy Implant - MIS Implants - and more	X		X	X	X	X		X	X	X	X	X	X	- CE	

**3Shape Headquarters**  
Europe, Middle East & Africa  
Holmens Kanal 7  
1060 Copenhagen, Denmark  
Tel: +45 7027 2620

**3Shape Asia**  
Room 906, Tower A of Eton Place  
No. 69, Dongfang Road  
200120 Shanghai, China  
Tel: +86 21 5835 2281

**3Shape Latin America**  
Carrera 13 # 82-91  
Oficina 401  
110221 Bogotá, Colombia  
Tel: +57 1691 9508

**3Shape North America**  
Somerset Hills Corporate Center  
10 Independence Boulevard, Suite 150  
Warren, New Jersey 07059, USA  
Tel: +1 (908) 867 0144

**More informat on at**  
www.3shape.com  
www.3shapedental.com  
info@3shape.com



Implant systems	Implant systems							Implant Solutions					Manu- facturing	Consumables (for local use)					Additional info				
	Providers of Implant solutions	Nobel Biocare	Biomet 3i	Straumann	Zimmer Dental	DENTSPLY - Astra Tech	DENTSPLY - Friudent	CAMLOG	Other systems	1- Piece (Titanium)	1- Piece (Zirconia)	2- Piece (Ti base + Zr abutment)	Implant Bars & Bridges	Digital Implant models	Providers milling center	Local milling by 3rd party	Blanks with pre-milled interface (for 1-piece abutment)	Titanium bases (for 2-piece abutments)	Intraoral Scan bodies	Digital model analogs	Multi-unit abutments	Cleared for use in US	Comments
GC Advanced Technologies (US) www.gc-a.com	reseller	👍	👍	👍	👍	👍			BioHorizons Sybron	X	X	X			X								
Glidewell (US) www.glidewell.com	reseller	👍	👍	👍	👍	👍	👍		Prisma & DentalCraft Keystone Dental Neoss	X	X	X	X	X			X	X	X		X		- CE pending
Green DenTech / Denrade (TW) www.denrade.com		👍	👍	👍	👍	👍	👍		Biohorizons Dentium Osstem and more...	X						X	X						- CE
Heraeus Kulzer (DE) www.heraeus.com	reseller	👍	👍	👍	👍	👍	👍		Thommen	X	X	X		X	X								
Implant Direct (US) www.implantdirect.com/custom-direct		👍		👍	👍				Biohorizons	X			X	X						X		X	- CE Health Canada
Intra-Lock International (US) www.intra-lock.com									INTRA-LOCK	X		X	X			X	X	X	X	X	X	X	- CE Health Canada
Ivoclar Wieland (DE) www.wieland-dental.de	reseller	👍	👍	👍	👍	👍	👍		Dentaurum Implants	X		X	X		X	X	X		X				
Lasak (CZ) www.lasak.com		👍		👍		👍			Lasak BlomIQ Lasak Implantent	X	X	X	X		X		X	X			X		- CE - ISO 13485 - ISO 9001
LaStruttura (IT) www.lastruttura.it	reseller	👍	👍	👍	👍	👍	👍		Sweden & Martina Magagn Prodent and many more	X	X		X	X				X	X				- CE 93/42
Medentika (DE) www.medentika.de		👍	👍	👍	👍	👍	👍		Medentika M-Implant	X		X	X		X	X	X	X			X		- CE
Medentis Medical (DE) www.medentis.de		👍		👍		👍	👍		ICX Templant	X		X	X	X			X	X					- CE
Medical Instinct (DE) www.medical-instinct.de									Medical Instinct	X	X	X	X		X	X	X	X			X		- CE
Médical Production (FR) www.medical-production.eu		👍	👍	👍	👍	👍	👍		Euroteknika	X		X	X	X	X	X	X	X			X		- CE0499 - ISO 9001 - ISO 13485
Megagen Implant (KR) www.imegagen.com									AnyBridge Internal AnyOne Internal MINI Internal Octa Level			X		X	X		X	X	X		X		- CE - KFDA
MIS Implants Technologies Ltd. (IL) - FDA cleared only www.mis-implants.com									MIS Implants			X	X		X	X	X	X			X	X	- CE
MIS Implants Technologies Ltd. (IL) - all libraries www.mis-implants.com									MIS Implants			X	X		X	X	X	X			X		- CE
Neodent (BR) www.neodent.com.br		👍		👍					Neodent Implant Systems	X	X	X	X		X			X				X	
Neoss (UK) - all libraries www.neoss.com		👍	👍	👍	👍	👍			Neoss	X		X	X	X	X	X	X	X	X	X	X	X	- CE
Nobel Biocare (CH) www.nobelbiocare.com		👍	👍	👍	👍	👍				X	X	X		X								X	- CE
NT Trading (DE) www.nt-trading.com		👍	👍	👍	👍	👍	👍		Thommen Sweden & Martina BEGO	X		X	X		X	X	X	X	X			X	- CE CMCAS GOST R TGC - Mexico Certification
Phibo (ES) www.phibo.com	reseller	👍	👍	👍	👍	👍	👍		BTI Sweden & Martina MIS Implants and many more	X	X	X	X	X				X	X				- CE - ISO 13485 - ISO 9001
Prowital (DE) www.prowital.de									Prowital			X			X		X	X	X				- CE - EN-ISO 9001-V472008 - EN-ISO-13485-2003AC2009 - Certificate Directive-93-42 EWG

**3Shape Headquarters**  
Europe, Middle East & Africa  
Holmens Kanal 7  
1060 Copenhagen, Denmark  
Tel: +45 7027 2620

**3Shape Asia**  
Room 906, Tower A of Eton Place  
No. 69, Dongfang Road  
200120 Shanghai, China  
Tel: +86 21 5835 2281

**3Shape Latin America**  
Carrera 13 # 82-91  
Oficina 401  
110221 Bogotá, Colombia  
Tel: +57 1691 9508

**3Shape North America**  
Somerset Hills Corporate Center  
10 Independence Boulevard, Suite 150  
Warren, New Jersey 07059, USA  
Tel: +1 (908) 867 0144

**More informat on at**  
www.3shape.com  
www.3shapedental.com  
info@3shape.com



Implant systems  Providers of Implant solutions	Implant systems							Implant Solutions					Manu- facturing		Consumables (for local use)						Additional info	
	Nobel Biocare	Biomet 3i	Straumann	Zimmer Dental	DENTSPLY - Astra Tech	DENTSPLY - Friudent	CAMLOG	Other systems	1- Piece (Titanium)	1- Piece (Zirconia)	2- Piece (Ti base + Zr abutment)	Implant Bars & Bridges	Digital Implant models	Providers milling center	Local milling by 3rd party	Blanks with pre-milled interface (for 1-piece abutment)	Titanium bases (for 2-piece abutments)	Intraoral Scan bodies	Digital model analogs	Multi-unit abutments	Cleared for use in US	Comments
Ritter Implants (DE) www.ritterimplants.com			👍	👍			- Alpha Bio Tec - MS Implants	X		X		X		X		X	X	X				- CE
Straumann (CH) - FDA cleared only www.straumann.com			👍					X	X	X			X								X	- CE
Straumann (CH) - all libraries www.straumann.com			👍					X	X	X			X	X	X	X	X	X	X			- CE
Sweden&Martina S.p.A (IT) www.sweden-martina.com							- Sweden&Martina	X	X	X	X	X	X			X	X	X				
Target3D (US) www.target3d.com	👍	👍	👍	👍	👍	👍	- Osstem - BioHorizons - MS-Implants - and more	X	X	X	X	X	X	X	X	X	X	X	X	X		- CE
T.F.I. System srl (IT) www.tfisystem.it							- Easy Grip			X	X			X		X	X	X	X	X	X	- CE
Thommen Medical (CH) www.thommenmedical.com							- Thommen			X	X					X	X				X	- CE
TRI Implants (CH) www.tri-implants.com							- TRI Implants			X						X	X	X				
Vulcan Custom Dental (US) www.vulcandental.com				👍			- BioHorizons	X		X	X	X	X	X	X	X	X	X	X	X	X	
Zfx (DE) www.zfx-dental.com	👍	👍	👍	👍	👍	👍	- Megagen - Thommen Medical - Microdent - Osstem - and more	X	X	X	X	X	X	X		X			X			- CE
Zimmer Dental (US) www.zimmerdental.com				👍				X	X				X								X	- CE

**Please note:**

- The above information is solely provided by third parties, and 3Shape assumes no responsibility for its accuracy or completeness.

DEPARTMENT OF HEALTH AND HUMAN SERVICES  
Food and Drug Administration**Indications for Use**

Form Approved: OMB No. 0910-0120

Expiration Date: January 31, 2017

See PRA Statement below.

510(k) Number (if known)

K151455

Device Name

3Shape Abutment Designer™ Software

## Indications for Use (Describe)

The 3Shape Abutment Designer Software is intended as an aid to the restoration of chewing function in partially or fully edentulous mandibles and maxillae. The 3Shape Abutment Designer Software is intended for use by a dental practitioner or dental laboratory staff for designing the patient specific component of a two-piece, one-piece, or hybrid dental implant abutment. The single or multi-unit abutment design is intended to be used by the manufacturer of an endosseous dental implant abutment to create the final device.

Type of Use (Select one or both, as applicable)

 Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

## **510(K) SUMMARY – Traditional 510(K)**

### **Submitter Information**

A Company Name: 3Shape A/S

B Company Address: Holmens Kanal 7  
DK-1060 Copenhagen K

C Company Phone: +45 7027 2620  
Company Fax: +45 7027 2621

D Contact Person: Hanne Nielsen  
Regulatory Affairs Manager

E Date Summary Prepared: June 28, 2016

### **Device Identification**

A Trade/proprietary Name: 3Shape Abutment Designer™  
Software

B Common Name: Abutment design software for dental  
laboratory

C Device Classification Name: Endosseous Dental Implant Abutment

C Regulation Number: 872.3630

C Classification: Class II

D Product Code: PNP

### **Predicate Device**

Sirona Dental CAD/CAM System (K100152).

### **Intended Use**

The 3Shape Abutment Designer Software is intended as an aid to the restoration of chewing function in partially or fully edentulous mandibles and maxillae. The 3Shape Abutment Designer Software is intended for use by a dental practitioner or dental laboratory staff for designing the patient specific component of a two-piece, one-piece or hybrid dental

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

implant abutment. The single or multi-unit abutment design is intended to be used by the manufacturer of an endosseous dental implant abutment to create the final device.

### **Device Description**

The 3Shape Abutment Designer™ Software receives scan data containing topographical characteristics of real teeth, dental impressions, or stone models. The software provides the user with the ability to create matching endosseous dental implant abutments using computer aided design.

The output of the device is a computer file containing the abutment(s) in digital form which can be used by manufacturers that hold an implant abutment 510(k) or are milling per the specific instructions provided by the holder of a 510(k) marketing clearance for a patient specific implant abutment. The 3Shape Abutment Designer Software output is restricted in the U.S. to be manufactured by a holder of an implant abutment 510(k) or milled per the specific instructions provided by the holder of a 510(k) marketing clearance for a patient specific implant abutment.

The 3Shape Abutment Designer™ Software can be run on properly configured "off-the-shelf" PC hardware running Microsoft® Windows and using a standard keyboard and mouse.

The 3Shape Abutment Designer™ Software requires the loading of implant libraries, which includes information such as implant type, maximum and minimum dimensional parameters for abutments, etc., created by separate abutment manufacturers and cleared by the FDA. In the US, the Implant Libraries are obtained via a 3Shape server after demonstration to 3Shape of the FDA clearance of the Implant Library.

### **Summary of the technological characteristics**

The 3Shape Abutment Designer™ is a software only device programmed in the Delphi programming language and has the following PC/laptop requirements and other technological characteristics as compared to the predicate device:



## 3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

Item	Submission Device Minimum Requirements	Submission Device Recommended	Predicate Device
OS	Windows 7 32-bit Professional*	Windows 7 64-bit Professional	Windows Vista 32-bit
RAM	4GB	8GB (16GB)	6GB
Video Card	512MB DirectX 10 (1GB DirectX 10)	1GB DirectX 11 (2GB DirectX 11)	512 MB DirectX 10 NVIDIA Quadro
Available HDD Space	250GB	500GB (1TB)	500GB
CPU	Intel Core i5 or equivalent	Intel Core i7 or equivalent	Intel Core i7
Monitor resolution	1440 x 900 pixels	1920 x 1080 pixels	Unknown
3D Mouse	None	3DConnexion SpaceMouse™	Unknown
Network	Internet connection		Unknown
Ports	USB 2.0 port for 3Shape desktop scanner		Serial Port PCI card for inEOS scanner
Mouse	Mouse with wheel button support		Unknown

## 3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

<b>Indications for Use</b>	The 3Shape Abutment Designer Software is intended as an aid in the restoration of chewing function in partially or fully edentulous mandibles and maxillae. The 3Shape Abutment Designer Software is intended for use by a dental practitioner or dental laboratory staff for designing the patient specific component of a two-piece, one-piece or hybrid dental implant abutment. The single or multi-unit abutment design is intended to be used by the manufacturer of an endosseous dental implant abutment to create the final device.	The Sirona Dental CAD/CAM System is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multi-unit cement retained restorations. The system consists of three major parts: TiBase, InCoris mesostructure, and CAD/CAM software. Specifically, the InCoris mesostructure and TiBase components make up a two-piece abutment which is used in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity. The InCoris mesostructure may also be used in conjunction with the CAMlog Titanium base CAD/CAM (types K2244.xxxx) (K083496) in the Camlog Implant System. The CAD/CAM software is intended to design and fabricate the InCoris mesostructure. The InCoris mesostructure and TiBase two-piece abutment is compatible with the following implant systems: Nobel Biocare Replace (K020646), Nobel Biocare Branemark (K022562), Friadent Xive (K013867), Biomet 3i Osseotite (K980549), Astra Tech Osseospeed (K091239), Zimmer Tapered Screw-Vent (K061410), and Straumann SynOcta (K061176).
<b>Software Output</b>	Digital encrypted or non-encrypted proprietary or .STL file including only the patient-specific abutment component for one-piece, two-piece, or hybrid abutment designs. The digital output does not include the abutment-to-implant connection interface.	.STL file of the ceramic mesostructure sent to Sirona Dental CAD/CAM System milling unit
<b>Physical - Output</b>	N/A – Submission device relies on separate regulatory clearance and manufacture of the abutment by a separate company	Two-piece TiBase abutment – pre-milled titanium base combined with ceramic mesostructure designed in Sirona CAD/CAM software.
<b>Milling Location</b>	Abutment Manufacturer or Dental laboratory per the 510(k) clearance of the dental abutment	Local milling of the ceramic abutment component.

The predicate Sirona device includes in the marketing clearance the two-piece TiBase Abutments (titanium bases and ceramic blocks) as a physical output as well as a validated milling unit and directions for assembly of the final dental abutment. The 3Shape Abutment Designer™ Software does not provide any physical parts that can come into contact with the patient. The 3Shape Abutment Designer Software instead provides only the digital design as an accessory to the physical dental abutment systems cleared by other manufacturers. The differences between the Indications for Use

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

Statement of the submission device and the predicate are related to the lack of physical output from the submission device and reliance on previous or subsequent FDA clearance of the physical abutment by a separate manufacturer. However, the difference does not change the intended use of the device.

### **Nonclinical Testing**

Software, hardware, and integration verification and validation testing was performed in accordance with the FDA Guidance Document "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices" (Issued on May 11, 2005) as well as the FDA Guidance Document "Off-The-Shelf Software Use in Medical Devices (Issued on September 9, 1999).

Prior to release, verification and validation testing of the 3Shape Abutment Designer™ Software has been completed using the approved acceptance criteria: Each user need has its own validation acceptance criteria; each specification has its own verification acceptance criteria; bug verification consists in ensuring issue is not reproducible; issues reported by beta partners must be reviewed and handled appropriately; beta partners must report functionality to be an improvement over previous version.

The validation suite includes validation of implemented mitigations related to device hazards identified in the risk management procedures.

### **Conclusion**

Based on a comparison of intended use, principle of operations, features and technical data, and the verification/validation test results, the 3Shape Abutment Designer™ Software is found to be substantially equivalent with the Predicate Device.

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

## **510(K) SUMMARY – Traditional 510(K)**

### **Submitter Information**

A Company Name: 3Shape A/S

B Company Address: Holmens Kanal 7  
DK-1060 Copenhagen K

C Company Phone: +45 7027 2620  
Company Fax: +45 7027 2621

D Contact Person: Hanne Nielsen  
Regulatory Affairs Manager

E Date Summary Prepared: June 28, 2016

### **Device Identification**

A Trade/proprietary Name: 3Shape Abutment Designer™  
Software

B Common Name: Abutment design software for dental  
laboratory

C Device Classification Name: Endosseous Dental Implant Abutment

C Regulation Number: 872.3630

C Classification: Class II

D Product Code: PNP

### **Predicate Device**

Sirona Dental CAD/CAM System (K100152).

### **Intended Use**

The 3Shape Abutment Designer Software is intended as an aid to the restoration of chewing function in partially or fully edentulous mandibles and maxillae. The 3Shape Abutment Designer Software is intended for use by a dental practitioner or dental laboratory staff for designing the patient specific component of a two-piece, one-piece or hybrid dental

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

implant abutment. The single or multi-unit abutment design is intended to be used by the manufacturer of an endosseous dental implant abutment to create the final device.

### **Device Description**

The 3Shape Abutment Designer™ Software receives scan data containing topographical characteristics of real teeth, dental impressions, or stone models. The software provides the user with the ability to create matching endosseous dental implant abutments using computer aided design.

The output of the device is a computer file containing the abutment(s) in digital form which can be used by manufacturers that hold an implant abutment 510(k) or are milling per the specific instructions provided by the holder of a 510(k) marketing clearance for a patient specific implant abutment. The 3Shape Abutment Designer Software output is restricted in the U.S. to be manufactured by a holder of an implant abutment 510(k) or milled per the specific instructions provided by the holder of a 510(k) marketing clearance for a patient specific implant abutment.

The 3Shape Abutment Designer™ Software can be run on properly configured "off-the-shelf" PC hardware running Microsoft® Windows and using a standard keyboard and mouse.

The 3Shape Abutment Designer™ Software requires the loading of implant libraries, which includes information such as implant type, maximum and minimum dimensional parameters for abutments, etc., created by separate abutment manufacturers and cleared by the FDA. In the US, the Implant Libraries are obtained via a 3Shape server after demonstration to 3Shape of the FDA clearance of the Implant Library.

### **Summary of the technological characteristics**

The 3Shape Abutment Designer™ is a software only device programmed in the Delphi programming language and has the following PC/laptop requirements and other technological characteristics as compared to the predicate device:

## 3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

Item	Submission Device Minimum Requirements	Submission Device Recommended	Predicate Device
OS	Windows 7 32-bit Professional*	Windows 7 64-bit Professional	Windows Vista 32-bit
RAM	4GB	8GB (16GB)	6GB
Video Card	512MB DirectX 10 (1GB DirectX 10)	1GB DirectX 11 (2GB DirectX 11)	512 MB DirectX 10 NVIDIA Quadro
Available HDD Space	250GB	500GB (1TB)	500GB
CPU	Intel Core i5 or equivalent	Intel Core i7 or equivalent	Intel Core i7
Monitor resolution	1440 x 900 pixels	1920 x 1080 pixels	Unknown
3D Mouse	None	3DConnexion SpaceMouse™	Unknown
Network	Internet connection		Unknown
Ports	USB 2.0 port for 3Shape desktop scanner		Serial Port PCI card for inEOS scanner
Mouse	Mouse with wheel button support		Unknown

## 3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

<b>Indications for Use</b>	The 3Shape Abutment Designer Software is intended as an aid in the restoration of chewing function in partially or fully edentulous mandibles and maxillae. The 3Shape Abutment Designer Software is intended for use by a dental practitioner or dental laboratory staff for designing the patient specific component of a two-piece, one-piece or hybrid dental implant abutment. The single or multi-unit abutment design is intended to be used by the manufacturer of an endosseous dental implant abutment to create the final device.	The Sirona Dental CAD/CAM System is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multi-unit cement retained restorations. The system consists of three major parts: TiBase, InCoris mesostructure, and CAD/CAM software. Specifically, the InCoris mesostructure and TiBase components make up a two-piece abutment which is used in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity. The InCoris mesostructure may also be used in conjunction with the CAMlog Titanium base CAD/CAM (types K2244.xxxx) (K083496) in the Camlog Implant System. The CAD/CAM software is intended to design and fabricate the InCoris mesostructure. The InCoris mesostructure and TiBase two-piece abutment is compatible with the following implant systems: Nobel Biocare Replace (K020646), Nobel Biocare Branemark (K022562), Friadent Xive (K013867), Biomet 3i Osseotite (K980549), Astra Tech Osseospeed (K091239), Zimmer Tapered Screw-Vent (K061410), and Straumann SynOcta (K061176).
<b>Software Output</b>	Digital encrypted or non-encrypted proprietary or .STL file including only the patient-specific abutment component for one-piece, two-piece, or hybrid abutment designs. The digital output does not include the abutment-to-implant connection interface.	.STL file of the ceramic mesostructure sent to Sirona Dental CAD/CAM System milling unit
<b>Physical - Output</b>	N/A – Submission device relies on separate regulatory clearance and manufacture of the abutment by a separate company	Two-piece TiBase abutment – pre-milled titanium base combined with ceramic mesostructure designed in Sirona CAD/CAM software.
<b>Milling Location</b>	Abutment Manufacturer or Dental laboratory per the 510(k) clearance of the dental abutment	Local milling of the ceramic abutment component.

The predicate Sirona device includes in the marketing clearance the two-piece TiBase Abutments (titanium bases and ceramic blocks) as a physical output as well as a validated milling unit and directions for assembly of the final dental abutment. The 3Shape Abutment Designer™ Software does not provide any physical parts that can come into contact with the patient. The 3Shape Abutment Designer Software instead provides only the digital design as an accessory to the physical dental abutment systems cleared by other manufacturers. The differences between the Indications for Use

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

Statement of the submission device and the predicate are related to the lack of physical output from the submission device and reliance on previous or subsequent FDA clearance of the physical abutment by a separate manufacturer. However, the difference does not change the intended use of the device.

### **Nonclinical Testing**

Software, hardware, and integration verification and validation testing was performed in accordance with the FDA Guidance Document "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices" (Issued on May 11, 2005) as well as the FDA Guidance Document "Off-The-Shelf Software Use in Medical Devices (Issued on September 9, 1999).

Prior to release, verification and validation testing of the 3Shape Abutment Designer™ Software has been completed using the approved acceptance criteria: Each user need has its own validation acceptance criteria; each specification has its own verification acceptance criteria; bug verification consists in ensuring issue is not reproducible; issues reported by beta partners must be reviewed and handled appropriately; beta partners must report functionality to be an improvement over previous version.

The validation suite includes validation of implemented mitigations related to device hazards identified in the risk management procedures.

### **Conclusion**

Based on a comparison of intended use, principle of operations, features and technical data, and the verification/validation test results, the 3Shape Abutment Designer™ Software is found to be substantially equivalent with the Predicate Device.



DEPARTMENT OF HEALTH AND HUMAN SERVICES  
Food and Drug Administration**Indications for Use**

Form Approved: OMB No. 0910-0120

Expiration Date: January 31, 2017

See PRA Statement below.

510(k) Number (if known)

K151455

Device Name

3Shape Abutment Designer™ Software

## Indications for Use (Describe)

The 3Shape Abutment Designer Software is intended as an aid to the restoration of chewing function in partially or fully edentulous mandibles and maxillae. The 3Shape Abutment Designer Software is intended for use by a dental practitioner or dental laboratory staff for designing the patient specific component of a two-piece, one-piece, or hybrid dental implant abutment. The single or multi-unit abutment design is intended to be used by the manufacturer of an endosseous dental implant abutment to create the final device.

Type of Use (Select one or both, as applicable)

 Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

DEPARTMENT OF HEALTH AND HUMAN SERVICES  
Food and Drug Administration

## Indications for Use

Form Approved: OMB No. 0910-0120

Expiration Date: January 31, 2017

See PRA Statement below.

510(k) Number (if known)

K151455

Device Name

3Shape Abutment Designer™ Software

Indications for Use (Describe)

The 3Shape Abutment Designer Software is intended as an aid to the restoration of chewing function in partially or fully edentulous mandibles and maxillae. The 3Shape Abutment Designer Software is intended for use by a dental practitioner or dental laboratory staff for designing the patient specific component of a two-piece, one-piece, or hybrid dental implant abutment. The single or multi-unit abutment design is intended to be used by the manufacturer of an endosseous dental implant abutment to create the final device.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRASStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

September 6, 2016

3Shape A/S  
Hanne Nielsen  
Regulatory Affairs Manager  
Holmens Kanal 7  
Copenhagen, 1060  
DENMARK

Re: K151455

Trade/Device Name: 3Shape Abutment Designer™ Software  
Regulation Number: 21 CFR 872.3630  
Regulation Name: Endosseous Dental Implant Abutment  
Regulatory Class: Class II  
Product Code: PNP  
Dated: January 8, 2016  
Received: January 11, 2016

Dear Hanne Nielsen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Page 2 - Hanne Nielsen

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

 Tina  
Kiang -S

Tina Kiang, Ph.D.  
Acting Director  
Division of Anesthesiology,  
General Hospital, Respiratory,  
Infection Control, and Dental Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES  
Food and Drug AdministrationForm Approved: OMB No. 0910-0120  
Expiration Date: January 31, 2017  
See PRA Statement below.**Indications for Use**

510(k) Number (if known)

K151455

Device Name

3Shape Abutment Designer™ Software

**Indications for Use (Describe)**

The 3Shape Abutment Designer Software is intended as an aid to the restoration of chewing function in partially or fully edentulous mandibles and maxillae. The 3Shape Abutment Designer Software is intended for use by a dental practitioner or dental laboratory staff for designing the patient specific component of a two-piece, one-piece, or hybrid dental implant abutment. The single or multi-unit abutment design is intended to be used by the manufacturer of an endosseous dental implant abutment to create the final device.

Type of Use (Select one or both, as applicable)

 Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

## **510(K) SUMMARY – Traditional 510(K)**

### **Submitter Information**

A Company Name: 3Shape A/S

B Company Address: Holmens Kanal 7  
DK-1060 Copenhagen K

C Company Phone: +45 7027 2620  
Company Fax: +45 7027 2621

D Contact Person: Hanne Nielsen  
Regulatory Affairs Manager

E Date Summary Prepared: June 28, 2016

### **Device Identification**

A Trade/proprietary Name: 3Shape Abutment Designer™  
Software

B Common Name: Abutment design software for dental  
laboratory

C Device Classification Name: Endosseous Dental Implant Abutment

C Regulation Number: 872.3630

C Classification: Class II

D Product Code: PNP

### **Predicate Device**

Sirona Dental CAD/CAM System (K100152).

### **Intended Use**

The 3Shape Abutment Designer Software is intended as an aid to the restoration of chewing function in partially or fully edentulous mandibles and maxillae. The 3Shape Abutment Designer Software is intended for use by a dental practitioner or dental laboratory staff for designing the patient specific component of a two-piece, one-piece or hybrid dental

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

implant abutment. The single or multi-unit abutment design is intended to be used by the manufacturer of an endosseous dental implant abutment to create the final device.

### **Device Description**

The 3Shape Abutment Designer™ Software receives scan data containing topographical characteristics of real teeth, dental impressions, or stone models. The software provides the user with the ability to create matching endosseous dental implant abutments using computer aided design.

The output of the device is a computer file containing the abutment(s) in digital form which can be used by manufacturers that hold an implant abutment 510(k) or are milling per the specific instructions provided by the holder of a 510(k) marketing clearance for a patient specific implant abutment. The 3Shape Abutment Designer Software output is restricted in the U.S. to be manufactured by a holder of an implant abutment 510(k) or milled per the specific instructions provided by the holder of a 510(k) marketing clearance for a patient specific implant abutment.

The 3Shape Abutment Designer™ Software can be run on properly configured "off-the-shelf" PC hardware running Microsoft® Windows and using a standard keyboard and mouse.

The 3Shape Abutment Designer™ Software requires the loading of implant libraries, which includes information such as implant type, maximum and minimum dimensional parameters for abutments, etc., created by separate abutment manufacturers and cleared by the FDA. In the US, the Implant Libraries are obtained via a 3Shape server after demonstration to 3Shape of the FDA clearance of the Implant Library.

### **Summary of the technological characteristics**

The 3Shape Abutment Designer™ is a software only device programmed in the Delphi programming language and has the following PC/laptop requirements and other technological characteristics as compared to the predicate device:

## 3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

Item	Submission Device Minimum Requirements	Submission Device Recommended	Predicate Device
OS	Windows 7 32-bit Professional*	Windows 7 64-bit Professional	Windows Vista 32-bit
RAM	4GB	8GB (16GB)	6GB
Video Card	512MB DirectX 10 (1GB DirectX 10)	1GB DirectX 11 (2GB DirectX 11)	512 MB DirectX 10 NVIDIA Quadro
Available HDD Space	250GB	500GB (1TB)	500GB
CPU	Intel Core i5 or equivalent	Intel Core i7 or equivalent	Intel Core i7
Monitor resolution	1440 x 900 pixels	1920 x 1080 pixels	Unknown
3D Mouse	None	3DConnexion SpaceMouse™	Unknown
Network	Internet connection		Unknown
Ports	USB 2.0 port for 3Shape desktop scanner		Serial Port PCI card for inEOS scanner
Mouse	Mouse with wheel button support		Unknown



## 3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

<b>Indications for Use</b>	The 3Shape Abutment Designer Software is intended as an aid in the restoration of chewing function in partially or fully edentulous mandibles and maxillae. The 3Shape Abutment Designer Software is intended for use by a dental practitioner or dental laboratory staff for designing the patient specific component of a two-piece, one-piece or hybrid dental implant abutment. The single or multi-unit abutment design is intended to be used by the manufacturer of an endosseous dental implant abutment to create the final device.	The Sirona Dental CAD/CAM System is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multi-unit cement retained restorations. The system consists of three major parts: TiBase, InCoris mesostructure, and CAD/CAM software. Specifically, the InCoris mesostructure and TiBase components make up a two-piece abutment which is used in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity. The InCoris mesostructure may also be used in conjunction with the CAMlog Titanium base CAD/CAM (types K2244.xxxx) (K083496) in the Camlog Implant System. The CAD/CAM software is intended to design and fabricate the InCoris mesostructure. The InCoris mesostructure and TiBase two-piece abutment is compatible with the following implant systems: Nobel Biocare Replace (K020646), Nobel Biocare Branemark (K022562), Friadent Xive (K013867), Biomet 3i Osseotite (K980549), Astra Tech Osseospeed (K091239), Zimmer Tapered Screw-Vent (K061410), and Straumann SynOcta (K061176).
<b>Software Output</b>	Digital encrypted or non-encrypted proprietary or .STL file including only the patient-specific abutment component for one-piece, two-piece, or hybrid abutment designs. The digital output does not include the abutment-to-implant connection interface.	.STL file of the ceramic mesostructure sent to Sirona Dental CAD/CAM System milling unit
<b>Physical - Output</b>	N/A – Submission device relies on separate regulatory clearance and manufacture of the abutment by a separate company	Two-piece TiBase abutment – pre-milled titanium base combined with ceramic mesostructure designed in Sirona CAD/CAM software.
<b>Milling Location</b>	Abutment Manufacturer or Dental laboratory per the 510(k) clearance of the dental abutment	Local milling of the ceramic abutment component.

The predicate Sirona device includes in the marketing clearance the two-piece TiBase Abutments (titanium bases and ceramic blocks) as a physical output as well as a validated milling unit and directions for assembly of the final dental abutment. The 3Shape Abutment Designer™ Software does not provide any physical parts that can come into contact with the patient. The 3Shape Abutment Designer Software instead provides only the digital design as an accessory to the physical dental abutment systems cleared by other manufacturers. The differences between the Indications for Use

3SHAPE ABUTMENT DESIGNER™ SOFTWARE 510(K) SUBMISSION – K151455

Statement of the submission device and the predicate are related to the lack of physical output from the submission device and reliance on previous or subsequent FDA clearance of the physical abutment by a separate manufacturer. However, the difference does not change the intended use of the device.

### **Nonclinical Testing**

Software, hardware, and integration verification and validation testing was performed in accordance with the FDA Guidance Document "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices" (Issued on May 11, 2005) as well as the FDA Guidance Document "Off-The-Shelf Software Use in Medical Devices (Issued on September 9, 1999).

Prior to release, verification and validation testing of the 3Shape Abutment Designer™ Software has been completed using the approved acceptance criteria: Each user need has its own validation acceptance criteria; each specification has its own verification acceptance criteria; bug verification consists in ensuring issue is not reproducible; issues reported by beta partners must be reviewed and handled appropriately; beta partners must report functionality to be an improvement over previous version.

The validation suite includes validation of implemented mitigations related to device hazards identified in the risk management procedures.

### **Conclusion**

Based on a comparison of intended use, principle of operations, features and technical data, and the verification/validation test results, the 3Shape Abutment Designer™ Software is found to be substantially equivalent with the Predicate Device.

## Substantial Equivalence Comparison

### 1. Predicates

The 3Shape Abutment Designer™ Software has the same intended uses and technical characteristics as the Sirona Dental CAD/CAM System (K100152) as listed in “Table 1: Predicate”

**Table 1: Predicate**

<b>Predicate</b>	<b>Manufacturer</b>	<b>510(k) number</b>	<b>Product code</b>
Dental CAD/CAM System	Sirona	K100152	NHA*

\* Endosseous dental implant abutments, 21CFR872.3630

## 2. Intended Use Comparison

### 2.1. 3Shape Abutment Designer™ Software

The Device's Intended Use, Intended users, and Intended Operational is stated in the "VOL\_001\_Administrative Documents" volume of this submission and reproduced here:

*"The 3Shape Abutment Designer™ Software is intended for designing endosseous dental implant abutments based on digital scans of 3D models.*

*Intended users are Dental Practitioners and Dental Laboratory staff.*

*Intended Operational Environment is Dental Laboratories."*

In the following sections, the similarities with the predicate are discussed.

### 2.2. Dental CAD/CAM System

The predicate's Intended Use can be extracted from the FDA 510(k) Premarket Notification Database:

*"The Sirona Dental CAD/CAM System is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations. The system consists of three major parts: TiBase, InCoris mesostructure, and CAD/CAM software. Specifically, the InCoris mesostructure and TiBase components make up a two-piece abutment which is used in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity.*

*The InCoris mesostructure may also be used in conjunction with the Camlog Titanium base CAD/CAM (types K2244.xxxx) (K083496) in the Canmlog Implant System. The CAD/CAM software is intended to design and fabricate the InCoris mesostructure. The InCoris mesostructure and TiBase two-piece abutment is compatible with the following implants systems:*

*[List of compatible implant systems omitted for brevity]*

Underlined segments indicate similarity with 3Shape Abutment Designer™ Software.

Note, the predicate is a CAD/CAM System bundled with physical two-piece abutments; the 3Shape Abutment Designer™ Software is CAD/CAM only.

### 3. Characteristics Comparison

A Characteristics Comparison can be seen in "Table 2: Substantial Equivalence Chart".

**Table 2: Substantial Equivalence Chart**

Feature name	3Shape Abutment Designer™ Software	Sirona Dental CAD/CAM System (K100152)	Detailed information
Graphical UI	Yes	Yes	See section 3.1
Windows OS platform	Yes	Yes	See section 3.2
Uses standard PC hardware	Yes	Yes	See section 3.3
Main software components	Order Management, Scan, Design, Production	Administration, Scan, Model, Design, Production	See section 3.4
Input data: Digitally imports topography of teeth by 3D Scan	Yes	Yes	See section 3.5
3D CAD design tools	Yes	Yes	See section 3.6
Custom abutment design	Yes	Yes	See intended use discussion in section 2.
Screw retained designs	Yes	Yes	See section 3.7
Implant Bridge design	Yes	Yes	See section 3.8
Patient safety measurements	Yes	Yes	See section 3.9
Export to milling	Yes	Yes	See section 3.10
Configurations and settings	Locked for end-users in the US	Proprietary information	See section 3.11
Technical documentation for 3 <sup>rd</sup> party software manufacturers	None	Proprietary information	N/A
Intended users	Dental practitioners and dental labs	Dental Technician or dentist	See section 3.12.
Output type	Computer file containing CAD model	Computer file containing CAD model	See section 3.10
Device submission includes pre-manufactured prosthetics*	No	Yes	See 510(k) summary of both devices.

\* Endosseous dental implant abutments as per 21CFR872.3630

**Note:** The comparison is made using the current publicly available labelling of Sirona Dental CAD/CAM System from the company's website (<http://manuals.sirona.com/en/digital-dentistry/>). The Predicate 510(k) summary lists the following components that is used for software comparison: CEREC MCXL, and inLab MCXL (labelling indicates rebranding as CEREC SW / CEREC SW Premium and inLab SW).

### 3.1. Graphical UI


Both devices utilizes a Graphical User Interface with a large 3D-based main window, a workflow progress bar, and various tools to manipulate the shape of the designed dental prosthetic / abutment.

Predicate device:



**Figure 1: Sirona Dental CAD/CAM System CEREC PREMIUM SW 4.4.x, -USA only-, 02.2016, Version: 121847**

3Shape Abutment Designer™:



**Figure 2: 3Shape Abutment Designer User interface as specified in requirements (See DS2015-1-RS0100 in the Additional Information 1 part of this submission).**

### 3.2. OS Platform

Both the predicate and the device run on the Windows OS platform.

#### Predicate Device:




Figure 3: Sirona Dental CAD/CAM System CEREC PREMIUM SW 4.4.x, -USA only-, 02.2016, Version: 121847

#### 3Shape Abutment Designer™:

3Shape Abutment Designer™ is compatible with Microsoft Windows 7 (64bit) /Microsoft Windows 8 (64bit) as specified in the *Architecture Design Chart* of the original submission.

### 3.3. Hardware Platform

Both devices runs on standard consumer PCs running Microsoft® Windows.

#### Predicate Device:

The Predicate Device fails to specify system requirements in the publically available labelling. However, the labelling repeatedly refers to the "PC", DVD and USB mediums, and the *Setup.exe* file (which is a format specific to PCs running Microsoft® Windows OS). The license is read from a USB stick attached to the PC.




Figure 4: Sirona Dental CAD/CAM System CEREC PREMIUM SW 4.4.x, -USA only-, 02.2016, Version: 121847





### 3Shape Abutment Designer™:

3Shape Abutment Designer™ is compatible with Microsoft® Windows 7 (64bit) / Microsoft® Windows 8 (64bit) as specified in the *Architecture Design Chart* of the original submission. The license is read from a USB dongle attached to a PC.

### 3.4. Main software components

Both the device and the predicate device are part of a larger system, which include Order Management, scan acquisition, General CAD tools, and Manufacturing interfaces. All these components are classified under product code NOF (510(k) exempt) with the exception of the Abutment design software.

The workflow and architecture of the two devices are almost identical as shown below. Note, the precise software architecture of the predicate device is proprietary information.

### Predicate Device:

The Software part of the Predicate comes in two flavours: CEREC™ targeted for Dentists / chair-side milling and InLab for Dental Laboratories. Publically available labelling for the two devices are essentially the same (see <http://manuals.sirona.com/en/digital-dentistry/>).

The predicate devices uses the following components of software represented by the workflow:

Substantial Equivalence Comparison.docx - Word



**Figure 5: Sirona Dental CAD/CAM System inLab SW 15.0, -USA only-, 08-2015, version: 120975**

### 3Shape Abutment Designer™:

As specified in the *Architecture Design Chart* of the original submission, the system, of which the device is a part of, consist of a number of components that match the following workflow:

- Dental Manager (Order Management / Administration)
- ScanIt (3D Scan acquisition, also includes the "Model" part of the Predicate Device)
- Dental Designer (CAD design)
- Manufacturing Interfaces (Production)



### 3.5. Input Data

Both devices rely on imports of digital topography of teeth produced by 3D Scanners. Both device manufacturers also market 3D Scanners.

Both devices use proprietary file formats, but can import the internationally recognized STL-format.

#### Predicate Device:

The Predicate Device labelling refers to the inEos line of 3D scanners for acquisition of the digital topography of teeth (refer to *Sirona Dental CAD/CAM System inLab SW 15.0, -USA only-, 08-2015, version: 120975* on [manufacturer website](#)).

The file formats used internally by the predicate is proprietary information and not available to 3Shape. However, the labelling states "*models or restorations can also be exported in \*.stl format for the further processing of this data in other software*".

#### 3Shape Abutment Designer™:

The device output is a 3D surface file in proprietary or STL-format.

### 3.6. 3D CAD design tools

Both the Device and the Predicate are part of a system with a long workflow process where at some point 3D CAD Design tools are used to design an abutment.

In the Predicate device, it happens in the DESIGN phase of the system. In 3Shape Abutment Designer™ this is done in Abutment Designer™ itself.

#### Predicate Device:



Figure 6: Example of 3D manipulation tools from Sirona Dental CAD/CAM System inLab SW 15.0, -USA only-, 08-2015, version: 120975





Figure 7: Example of 3D manipulation tools from Sirona Dental CAD/CAM System inLab SW 15.0, -USA only-, 08-2015, version: 120975

**3Shape Abutment Designer™:**

**►Step 4: Design the Abutment**

- Use the available control points and settings to design the abutment top cap and emergence profile.
- Right-click on the abutment and select **Add profile** to add additional vertical profiles with points for a better manipulation of the whole abutment. Right-click the profile to remove it.
- Right-click on the emergence profile and select **Add curvature point** to place extra horizontal profiles with points for a better manipulation of the emergence profile. Right-click the point to remove it.



**Robotic Abutment Settings:**

**Hint!** If necessary, you can change the abutment settings (e.g. an Implant Kit) in the Order form by clicking the corresponding button in the Workflow bar.

7

**Figure 8: Example of 3D manipulation tools from DS-2.15.2.0-A-EN\_Abutment\_Designer\_Excerpt\_Version included in the original submission**

**Parametric Customised Abutment**

Settings

Type: Robotic Abutment

Draft angle: 3.0 deg

Snap gingiva:

Visual options:

Advanced

Shoulder radius: 0.7 mm

Top fillet radius: 0.3 mm

Margin fillet radius: 0.00 mm

Grid size: 5.4 mm

Vertical offset: 0.00 mm

Neck wall thickness: 0.0 mm

Neck fillet: 0.0 mm

**Draft angle** - the side angle of the upper part of a top cap.

**Lock** - when marked, locks the draft angle during manual editing of the selected abutment.

**Snap gingiva** - automatically places abutment margin line on gingiva.

**Modify margin shape** - opens a window (see image) to select a suitable margin shape with the purpose to initially match the gingiva before manual editing of the abutment margin line.

**Save** - saves modified values as a default for the current and symmetrical abutments.

**Reset** - resets modified unsaved values to default.

**Visual options:**

- Shows/Hides angle graphics (insertion direction and implant direction).
- Shows/Hides abutment screw hole

**Advanced:**

- The **Edit Implant Models** button, in **Advanced settings**, lets you rotate the abutment base and position the axis of the screw hole at an angle with the cursor.
- The **Add bevel on shoulder** button lets you select the **Bevel on shoulder** option to create an additional surface above the emergence profile with the **Angle** and **Length** parameters set in relation to the top cap direction.
- **Shoulder radius** - the rounding radius of the lower part of a top cap.
- **Top fillet radius** - the rounding radius of the upper part of a top cap.
- **Margin fillet radius** - the rounding radius of abutment near the margin line.

8

**Figure 9: Example of 3D manipulation tools from DS-2.15.2.0-A-EN\_Abutment\_Designer\_Excerpt\_Version included in the original submission**



### 3.7. Screw Retained Designs

Both the Predicate Device and the Device offers abutments that are fixed to the implant with a screw. There is some subtle differences in how the design tools work, but the basic concept is the same.

#### Predicate Device:

The 510(k) summary states *"The two piece abutment is mounted onto the implant and fixed with a screw."*

#### 3Shape Abutment Designer™:

Screw retained design is part of the requirements as described in the SRS (see DS2015-1-RS0121 of the original submission).

### 3.8. Implant Bridge Design

Both devices offer the user the possibility of designing a Bridge restoration that can be mounted on top of implant abutments.

#### Predicate Device:

See section 8.7: *"Anatomical or reduced directly screwed-on bridges"* in the publically available labelling: *"Sirona Dental CAD/CAM System inLab SW 15.0, -USA only-, 08-2015, version: 120975"*

#### 3Shape Abutment Designer™:

See requirements in the SRS (DS2015-1-RS0127) and section 1.18: *"Implant Bars and Bridges"* of *"DS-2.15.2.0-A-EN\_Abutment\_Designer\_Excerpt\_Version"* included in the original submission.

### 3.9. Patient Safety Measures

The safety measures of this device includes and expands on those of the predicate device's. Where the predicate device only detects the angle of the designed abutment, 3Shape Abutment Designer™ checks a number of safety limitations as described below.


Where the predicate device only warns the user, 3Shape Abutment Designer™ forces a hard stop in the design process.

#### Predicate Device:

When designing abutments on the Predicate Device, the user receives a warning when the abutment exceeds the specified safety limits, but the user is allowed to continue designing.

From the official labelling on the Predicate Device, the device only has an Abutment Angle warning.

### 9.1.9 Define restoration axis



The angle between the implant axis (dotted red) and the restoration axis (yellow arrow) may be no higher than 20° (A).

If an angle of more than 20° is selected between the implant axis (dotted red) and the restoration axis (yellow arrow), the user will be warned by a change of color to red (B).

#### CAUTION

Exceeding specified safety limits of your device results in the construction of a misbranded device which may lead to premature abutment fracture. In this case the patient must be informed that he is going to receive a device that is beyond the labeled specifications. Please refer to the corresponding manual of the used material for the cleared angle/wall thickness limitations.

- A caution will pop up:  
"Warning! The angle you have set between the implant direction and abutment direction leads to a construction that is beyond the cleared safety limit. If you do not adjust this situation you are producing a device beyond the labeled specifications. Please consider the user manual for required actions".

**Figure 10 Sirona Dental CAD/CAM System CEREC PREMIUM SW 4.4.x, -USA only-, 02.2016, Version: 121847**

#### 3Shape Abutment Designer™:


When designing an abutment using Abutment Designer™, the user experiences both a warning and a 'hard stop', meaning that the user cannot continue designing if the abutment design exceeds the safety limits (DS2015-1-RS0121 and DS2015-1-RS0127 from the original submission).

3Shape Abutment Designer™ also holds some safety measures in the configuration and settings (see Section 3.11).

#### 3.10. Export to milling (output type)

Both devices are capable of exporting the digital file containing the completed abutment to a milling machine (**Note**, the Predicate Device manufacturer also markets milling machinery, 3Shape does not) either by means of a proprietary format or the STL-format.

#### Predicate Device:



**Figure 11: Sirona Dental CAD/CAM System inLab SW 15.0, -USA only-, 08-2015, version: 120975**

### 3Shape Abutment Designer™:

See "Device Description" in the submitted 510(k) summary:

*"The output of the device is a computer file containing the abutment(s) in digital form which can be used by 3rd party providers of milling machines and/or rapid prototyping machines to manufacture physical abutments."*

### 3.11. Configuration and settings

The 3Shape Device does not allow changing of Implant settings by the end-user for Implant Libraries cleared for use in the United States.


The Predicate device's handling of configurations is not discernible from the publically available labelling. However, labelling contains warnings that may be relevant.

### Predicate Device:



From the official labelling of the Predicate Device the user can change the configuration and is advised only to consult the manual of the individual implant system.

### 9.1.10 Adjusting parameters



- ✓ The step "*Restoration Parameters*" is active.
- 1. Adjust the parameters where necessary.
- 2. Confirm the changes with "*Ok*".
  - ↳ The initial suggestion is then calculated automatically.
- 3. If you have not altered any parameters, change to the "*Calculate Restoration*" step.
  - ↳ An initial suggestion is calculated.

**Figure 12 Sirona Dental CAD/CAM System CEREC PREMIUM SW 4.4.x, -USA only-, 02.2016, Version: 121847**

#### 3Shape Abutment Designer™:

3Shape Abutment Designer has a safety measure such that safety limits cannot be altered by the user, but is set by the implant provider (See DS2015-1-SS0012 and DS2015-1-SS0044).

### 3.12. Intended users

The intended users differ only in wording, "*dental technician or dentist*" vs. "*dental practitioners and dental laboratory staff*".

#### Predicate Device:

## 2 General data

Please read this document completely and follow the instructions exactly. You should always keep it within reach.

Original language of the present document: German

### 2.1 Certification

#### CE mark

This product bears the CE mark in accordance with the provisions of the Council Directive 93/42/EEC of June 14, 1993 concerning medical devices (MDD).



### 2.2 General safety information

#### Only use original software

Only use original software or software which has been released by Sirona. To produce restorations and equipment, manipulated or non-released software components must not be used.

Software and software components must not be installed using incorrect data.

Please check that each installed component has been granted approval in its country. Contact your dealer for more information.

#### Restoration to be checked by trained personnel

Each restoration which is performed with this software must be checked for suitability by a trained person (e.g. dental technician or dentist).

Figure 13: Sirona Dental CAD/CAM System inLab SW 15.0, -USA only-, 08-2015, version: 120975

### 3Shape Abutment Designer™:

See the Intended Use Comparison in section 2.



#### **4. Conclusion**

The 3Shape Abutment Designer™ Software and the predicate only deviate significantly in the cases where the predicate is bundled with a physical dental implant abutment. 3Shape Abutment Designer™ Software does not provide any physical parts that can come into contact with the patient.

The differences between the 3Shape Abutment Designer™ Software and the predicate do not raise additional concerns with respect to the safety and effectiveness of the 3Shape Abutment Designer™ Software.

Based on the information presented, we conclude Substantial Equivalence between the predicate and the 3Shape Abutment Designer™ Software.

# Standard Operating Procedure – SW-SOP-0011

TITLE: HANDLING OF FDA CLEARED IMPLANT LIBRARIES IN 3SHAPE SOFTWARE

VERSION: 2.0

APPROVAL STATUS: APPROVED

LAST CHANGED: 09-05-2016 07:59:00

LAST CHANGED BY: 3SHAPE\THOMAS CLEMEN

APPROVAL DATE: 09-05-2016 07:59:00

APPROVED BY: 3SHAPE\THOMAS CLEMEN

## 1 Introduction

3Shape Software Products are locked, so only implant libraries that have been cleared by the FDA can be used by users located in the United States.

An implant library is a set of digital versions of one or more implant systems. Clearing the implant library for use in the United States, will make all the implant systems in said library usable in the United States.

An implant library is created by an Implant Library Provider. The Implant Library Provider must create an implant library, such that it contains only implant systems and components that are cleared by the FDA. It is allowed for an Implant Provider to include implant systems originating from another provider, as long as the Implant System and all its components are cleared by the FDA.

This document describes the process of including an FDA cleared implant library, so it can be used in 3Shape Software Products.

## 2 Scope

This document applies to the Databank Engineers and Software Project Managers.

## 3 Description

Upon receiving a request from an Implant Library Provider to enable a specific implant library for use in the US, 3Shape must use the following procedure:



### 3.1 Create DMS Record

Each implant library, cleared for use in the US, must hold a corresponding DMS (Document Management System) record.

The DMS record must hold information about

- The name of the Implant Library Provider
- The regions in which this library should be cleared for use.

Use the following naming convention when creating the DMS record:

"Implant Library Clearance - <name of digital signature>.docx"

### **3.2 Review documentation from Implant Library Provider**

When 3shape receives a request from an Implant Library Provider about enabling an implant library for use in the US, it must contain:

1. A signed confirmation letter from the library provider stating that all components in the implant library are cleared by the FDA
2. FDA 510(k) clearance numbers(s) for the implant library.

If 3Shape holds an official record (such as a User Manual) listing the content of the implant libraries, then the Databank Engineer must ensure that the record corresponds to cleared implant library.

The documentation received from the Implant Library Providers must be stored in the DMS record.

### **3.3 Create Digital Signature**

The Databank Engineer must either update the existing digital signature used by that Implant Library Provider or create a new digital signature that the implant provider can use for encrypting the implant libraries.

The Databank Engineer must ensure setup of the digital signature so the regional setting on the digital signature includes USA.

The name of the digital signature must be stored in the DMS Record.

### **Step 4: Verify the library in the requested Software Product**

A verification of the encryption of the implant library must be performed and included into the DMS record.

### **Step 5: Approve the DMS Record**

The DMS record must be approved by a 3Shape Software Project Manager.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

September 6, 2016

3Shape A/S  
Hanne Nielsen  
Regulatory Affairs Manager  
Holmens Kanal 7  
Copenhagen, 1060  
DENMARK

Re: K151455

Trade/Device Name: 3Shape Abutment Designer™ Software  
Regulation Number: 21 CFR 872.3630  
Regulation Name: Endosseous Dental Implant Abutment  
Regulatory Class: Class II  
Product Code: PNP  
Dated: January 8, 2016  
Received: January 11, 2016

Dear Hanne Nielsen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Page 2 - Hanne Nielsen

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

 Tina  
Kiang -S

Tina Kiang, Ph.D.  
Acting Director  
Division of Anesthesiology,  
General Hospital, Respiratory,  
Infection Control, and Dental Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure



Contains Nonbinding Recommendations

Print Form

# Acceptance Checklist for Traditional 510(k)s

**(Should be completed within 15 days of DCC receipt)**  
The following information is not intended to serve as a comprehensive review.

510(k) #: K151455 Date Received by DCC: June 1, 2015  
Lead Reviewer: Andrew I. Steen - Trade Name: 3Shape Abutment Designer Software  
Branch: DEDB Division: DAGRID Center/Office: CDRH/ODE

Note: If an element is left blank on the checklist, it does not mean the checklist is incomplete. It means the reviewer did not assess the element during RTA and the element will be assessed during the substantive review.

Preliminary Questions		
Answers in the shaded blocks indicate consultations with Center advisor is needed	Yes	No
<p><b>1) Is the product a device (per section 201(h) of the FD&amp;C Act) or a combination product (per 21 CFR 3.2(e)) with a device constituent part subject to review in a 510(k)?</b></p> <p>If it appears not to be a device (per section 201(h) of the FD&amp;C Act) or such a combination product, or you are unsure, consult with the CDRH Jurisdictional Officer or the CBER Office Jurisdiction Liaison to determine the appropriate action, and inform division management. <i>Provide a summary of the Jurisdictional Officer's/Liaison's determination.</i> If the product does not appear to be a device or such a combination product, mark "No."</p>	X	
Comments?		
<p><b>2. Is the application with the appropriate Center?</b></p> <p>If the product is a device or a combination product with a device constituent part, is it subject to review by the Center in which the submission was received? If you believe the application is not with the appropriate Center or you are unsure, consult with the CDRH Jurisdictional Officer or CBER Office Jurisdiction Liaison to determine the appropriate action and inform your division management. <i>Provide a summary of the Jurisdictional Officer's/Liaison's determination.</i> If application should not be reviewed by your Center mark "No."</p>	X	
Comments?		
<p><b>3) If a Request for Designation was submitted for the device or combination product with a device constituent part and assigned to your center, identify the RFD # and confirm the following:</b></p> <p>a) Is the device or combination product the same (e.g., design, formulation) as that presented in the RFD submission?</p> <p>b) Are the indications for use for the device or combination product identified in the 510(k) the same as those identified in the RFD submission ?</p> <p>If you believe the product or the indications presented in the 510(k) have changed from the RFD, or you are unsure, consult with the CDRH Jurisdictional Officer or appropriate CBER Jurisdiction Liaison to determine the appropriate action and inform your division management. <i>Provide summary of Jurisdictional Officer's/Liaison's determination.</i> If the answer to either question is no, mark "No." If there was no RFD, skip this question.</p>		
Comments?		
<p><b>4) Is this device type eligible for a 510(k) submission?</b></p> <p>If a 510(k) does not appear to be appropriate (e.g., Class III type and PMA required, or Class I or II type and 510(k)-exempt), you should consult with the CDRH 510(k) Program Director or appropriate CBER staff during the acceptance review. If 510(k) is not the appropriate regulatory submission, mark "No."</p>	X	
Comments?		



<p>5) Is there a pending PMA for the same device with the same indications for use?</p>	<p>If yes, consult division management and the CDRH 510(k) Program Director or appropriate CBER staff to determine the appropriate action.</p>	<p>Comments?</p>
<p>6) If clinical studies have been submitted, is the submitter the subject of an Application Integrity Policy (AIP)?</p>	<p>If yes, consult with the CDRH Office of Compliance/Division of Bioresearch Monitoring (OC/DBM - BIMO) or CBER Office of Compliance and Biologics Quality/Division of Inspections and Surveillance/Bioresearch Monitoring Branch (OCBQ/DIS/BMB) to determine the appropriate action. Check on web at <a href="http://www.fda.gov/ICECI/EnforcementActions/ApplicationIntegrityPolicy/ucm134453.htm">http://www.fda.gov/ICECI/EnforcementActions/ApplicationIntegrityPolicy/ucm134453.htm</a></p>	<p>Comments?</p>

If the answer to 1 or 2 appears to be "No," then stop review of the 510(k) and issue the "Original Jurisdictional Product" letter.  
 If the answer to 3a or 3b appears to be "No," then stop the review and contact the CDRH Jurisdictional Officer or CBER Office of Jurisdiction Liaison.  
 If the answer to 4 is "No," the lead reviewer should consult division management and other Center resources to determine the appropriate action.  
 If the answer to 5 is "Yes," then stop review of the 510(k), contact the CDRH 510(k) Staff and PMA Staff, or appropriate CBER staff.  
 If the answer to 6 is "Yes," then contact CDRH/OC/DBM-BIMO or CBER/OCBQ/DIS/BMB, provide a summary of the discussion with the BIMO Staff, and indicate BIMO's recommendation/action.

**Organizational Elements**

*Failure to include these items alone generally should not result in an RTA designation.*

	Yes	No
1) Submission contains a Table of Contents	X	
2) Each section is labeled (e.g., headings or tabs designating Device Description section, Labeling section, etc.)	X	
3) All pages of the submission are numbered.	X	
4) Type of 510(k) is identified (i.e., traditional, abbreviated, or special)	X	
Comments?		

**Elements of a Complete Submission (RTA Items)****(21 CFR 807.87 unless otherwise indicated)**

Submission should be designated RTA if not addressed.

**Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.**

- Any "No" answer will result in a "Refuse to Accept" decision. - Each element on the checklist should be addressed within the submission. An applicant may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criteria is considered Present (Yes). An assessment of the rationale will be considered during the review of the submission.	Yes	No	N/A	Comment
---	-----	----	-----	---------

**A. Administrative**

1) All content used to support the submission is written in English (including translations of test reports, literature articles, etc.)	X			
2) Submission identifies the following (such as in CDRH Premarket Review Submission Cover Sheet (Form 3514) or 510(k) cover letter):	X			
a) Device trade name or proprietary name	X			
b) Device common name	X			
c) Device class and panel or Classification regulation or Statement that device has not been classified with rationale for that conclusion	X			
3) Submission contains Indications for Use Statement with Rx and/or OTC designation (see also 21 CFR 801.109).	X			
4) Submission contains 510(k) Summary or 510(k) Statement	X			
a) Summary contains all elements per 21 CFR 807.92 (See also 510(k) Summary Checklist)			X	
b) Statement contains all elements per 21 CFR 807.93			X	
5) Submission contains Truthful and Accuracy Statement per 21 CFR 807.87(k) See recommended format.	X			
6) Submission contains Class III Summary and Certification. See recommended content.			X	
7) Submission contains clinical data			X	
8) If submission references use of a national or international standard as part of demonstration of substantial equivalence, submission contains Standards Data Report for 510(k)s (Form 3654) or includes detailed information about how and the extent to which the standard has been followed.	X			
9) The submission identifies prior submissions for the same device for which FDA provided feedback related to the data or information needed to support substantial equivalence (e.g., submission numbers for Pre-Submission, IDE, prior not substantially equivalent (NSE) determination, prior 510(k) that was deleted or withdrawn) or states that there were no prior submissions for the subject device.	X			
a) If there were prior submissions, the submitter has identified where in the current submission any issues related to a determination of substantial equivalence outlined in prior communications are addressed. For additional information regarding the Pre-Submission process, please refer to the Draft Guidance "Medical Devices: The Pre-Submission Program and Meetings with FDA Staff." Once finalized, this guidance will represent the Agency's current thinking on this topic.	X			

**B. Device Description**

10)				
-----	--	--	--	--

**Elements of a Complete Submission (RTA Items)****(21 CFR 807.87 unless otherwise indicated)**

Submission should be designated RTA if not addressed.

**Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.**

- Any "No" answer will result in a "Refuse to Accept" decision. - Each element on the checklist should be addressed within the submission. An applicant may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criteria is considered Present (Yes). An assessment of the rationale will be considered during the review of the submission.	Yes	No	N/A	Comment
a) If there are requirements regarding the device description, such as special controls, in a device-specific regulation that are applicable to the device, the submission includes device description information to establish that the submitter has followed the device-specific requirement.			X	
b) If there is a device-specific guidance, other than a special controls guidance document, applicable to the device, the submission includes device description information to establish that the submitter has addressed the recommendations or otherwise has met the applicable statutory or regulatory criteria through an alternative approach.	X			
11) Descriptive information is present and consistent within the submission (e.g., the device description section is consistent with the device description in the labeling), including:				
a) A description of the principle of operation and mechanism of action for achieving the intended effect.	X			
b) A description of proposed conditions of use, such as surgical technique for implants; anatomical location of use; user interface; how the device interacts with other devices; and/or how the device interacts with the patient.	X			
c) A list and description of each device for which clearance is requested.	X			
12) Submission contains representative engineering drawing(s), schematics, illustrations and/or figures of the device that are clear, legible, labeled, and include dimensions.			X	
13) If device is intended to be marketed with multiple components, accessories, and/or as part of a system				
a) Submission includes a list of all components and accessories to be marketed with the subject device.	X			
b) Submission includes a description (as detailed in item 11(a) and (b) and 12 above) of each component or accessory.			X	
c) A 510(k) number is provided for each component or accessory that received a prior 510(k) clearance.			X	
<b>C. Substantial Equivalence Discussion</b>				
14) Submitter has identified a predicate device.	X			
a) Predicate's 510(k) number, trade name, and model number (if applicable) provided. For predicates that are preamendments devices, information is provided to document preamendments status. <i>Information regarding documenting preamendment status is available online.</i>	X			
b) The identified predicate(s) is consistent throughout the submission (i.e., the predicate(s) identified in the Substantial Equivalence section is the same as that listed in the 510(k) Summary (if applicable) and that used in comparative performance testing.	X			
15) Submission includes a comparison of the following for the predicate(s) and subject device				
a) Indications for Use	X			
b) Technology, including features, materials, and principles of operation	X			

**Elements of a Complete Submission (RTA Items)****(21 CFR 807.87 unless otherwise indicated)**

Submission should be designated RTA if not addressed.

**Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.**

- Any "No" answer will result in a "Refuse to Accept" decision. - Each element on the checklist should be addressed within the submission. An applicant may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criteria is considered Present (Yes). An assessment of the rationale will be considered during the review of the submission.	Yes	No	N/A	Comment
16) Submission includes an analysis of why any differences between the subject device and predicate(s) do not render the device NSE (e.g., does not constitute a new intended use; and any differences in technological characteristics are accompanied by information that demonstrates the device is as safe and effective as the predicate and do not raise different questions of safety and effectiveness than the predicate ), affect safety or effectiveness, or raise different questions of safety and effectiveness (see section 513(i)(1)(A) of the FD&C Act and <u>21 CFR 807.87(f)</u> )		X		X
Comments? Reviewer Discretion - the company states the differences do not raise additional concerns but does not explain how or why. This issue will be addressed during the review stage.				

**D. Proposed Labeling (see also 21 CFR part 801)**If *in vitro* diagnostic (IVD) device, criteria 17 & 19 may be omitted.

17) Submission includes proposed package labels and labeling (e.g., instructions for use, package insert, operator's manual) that include a description of the device, its intended use, and the directions for use.	X			X
a) Indications for use are stated in labeling and are identical to Indications for Use form and 510(k) Summary (if 510(k) Summary provided).		X		
b) Submission includes directions for use that - include statements of all conditions, purposes or uses for which the device is intended (e.g., hazards, warnings, precautions, contraindications) AND - includes directions for layperson (see <u>21 CFR 801.5</u> ) OR submission states that device qualifies for exemption per <u>21 CFR 801 Subpart D</u>		X		
Comments? Reviewer Discretion. This issue will be addressed during the review				
18) If indicated for prescription use, labeling includes the prescription use statement (see <u>21 CFR 801.109(b)(1)</u> ) or "Rx only" symbol [See also <u>Alternative to Certain Prescription Device Labeling Requirements</u> ]	X			
19) General labeling provisions				
a) Labeling includes name and place of business of the manufacturer, packer, or distributor (21 CFR 801.1).	X			
b) Labeling includes device common or usual name. (21 CFR 801.61)			X	
20)				
a) If there are requirements regarding labeling, such as special controls, in a device-specific regulation that are applicable to the device, the submission includes labeling to establish that the submitter has followed the device-specific requirement.			X	
b) If there is a device-specific guidance, other than a special controls guidance document, applicable to the device, the submission includes labeling to establish that the submitter has addressed the recommendations or otherwise has met the applicable statutory or regulatory criteria through an alternative approach.			X	

**Elements of a Complete Submission (RTA Items)****(21 CFR 807.87 unless otherwise indicated)**

Submission should be designated RTA if not addressed.

**Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.**

- Any "No" answer will result in a "Refuse to Accept" decision.

- Each element on the checklist should be addressed within the submission. An applicant may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criteria is considered Present (Yes). An assessment of the rationale will be considered during the review of the submission.

	Yes	No	N/A	Comment
c) If there is a special controls document applicable to the device, the submission includes labeling to establish that the submitter has complied with the particular mitigation measures set forth in the special controls document or uses alternative mitigation measures but provides a rationale to demonstrate that those alternative measures identified by the firm will provide at least an equivalent assurance of safety and effectiveness.	X			
21) If the device is an in vitro diagnostic device, provided labeling includes all applicable information required per <u>21 CFR 809.10</u> .			X	

**E. Sterilization**

If IVD device and sterilization is not applicable, select "N/A" and criteria below will be omitted from checklist.

Submission states that the device and/or accessories are: (one of the below must be checked)

provided sterile

provided non-sterile but sterilized by the end user

X non-sterile when used

Information regarding the sterility status of the device is not provided.

This information will determine whether and what type of additional information may be necessary for a substantial equivalence determination.

22) Assessment of the need for sterilization information

	Yes	No	N/A	Comment
a) Identification of device, and/or accessories, and/or components that are provided sterile.		X		
b) Identification of device, and/or accessories, and/or components that are end user sterilized.			X	
c) Identification of device, and/or accessories, and/or components that are reusable and cleaning /disinfection instructions are provided.			X	

25)

	Yes	No	N/A	Comment
a) If there are requirements regarding sterility, such as special controls, in a device-specific regulation that are applicable to the device, the submission includes sterility information to establish that the submitter has followed the device-specific requirement.			X	
b) If there is a device-specific guidance, other than a special controls guidance document, applicable to the device, the submission includes sterility information to establish that the submitter has addressed the recommendations or otherwise has met the applicable statutory or regulatory criteria through an alternative approach.			X	
c) If there is a special controls document applicable to the device, the submission includes sterility information to establish that the submitter has complied with the particular mitigation measures set forth in the special controls document or uses alternative mitigation measures but provides a rationale to demonstrate that those alternative measures identified by the firm will provide at least an equivalent assurance of safety and effectiveness.	X			

**Elements of a Complete Submission (RTA Items)****(21 CFR 807.87 unless otherwise indicated)**

Submission should be designated RTA if not addressed.

**Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.**

- Any "No" answer will result in a "Refuse to Accept" decision.  
 - Each element on the checklist should be addressed within the submission. An applicant may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criteria is considered Present (Yes). An assessment of the rationale will be considered during the review of the submission.

Yes

No

N/A

Comment

**F. Shelf Life**

26) Proposed shelf life/expiration date stated

X

27) For sterile device, submission includes summary of methods used to establish that device will remain sterile through the proposed shelf life or a rationale for why testing to establish shelf life is not applicable.

X

28) Submission includes summary of methods used to establish that device performance is not adversely affected by aging or includes a rationale for why the storage conditions are not expected to affect device safety or effectiveness.

X

**G. Biocompatibility**

If IVD device, select "N/A" and the below criteria will be omitted from checklist.

Submission states that there: (one of the below must be checked)

are direct or indirect (e.g., through fluid infusion) patient-contacting components.

X

are no direct or indirect (e.g., through fluid infusion) patient-contacting components.

Information regarding the patient contact status of the device is not provided.

This information will determine whether and what type of additional information may be necessary for a substantial equivalence determination.

**H. Software**

Submission states that the device: (one of the below must be checked)

X

does contain software/firmware.

does not contain software/firmware.

Information regarding whether the device contains software is not provided.

This information will determine whether and what type of additional information may be necessary for a substantial equivalence determination.

32) Submission includes a statement of software level of concern and rationale for the software level of concern.

X

33) All applicable software documentation provided based on level of concern identified by the submitter, as described in [Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices](#), or the submitter has provided an alternative approach with a rationale.

X

**I. EMC and Electrical Safety**

Submission states that the device: (one of the below must be checked)

does require EMC and Electrical Safety evaluation.

X

does not require EMC and Electrical Safety evaluation.

**Elements of a Complete Submission (RTA Items)**

**(21 CFR 807.87 unless otherwise indicated)**

Submission should be designated RTA if not addressed.

Check "Yes" if item is present, "N/A" if it is not needed and "No" if it is not included but needed.

- Any "No" answer will result in a "Refuse to Accept" decision. - Each element on the checklist should be addressed within the submission. An applicant may provide a rationale for omission for any criteria that are deemed not applicable. If a rationale is provided, the criteria is considered Present (Yes). An assessment of the rationale will be considered during the review of the submission.	Yes	No	N/A	Comment
---	-----	----	-----	---------

Information regarding whether the device requires EMC and Electrical Safety evaluation is not provided.

This information will determine whether and what type of additional information may be necessary for a substantial equivalence determination.

**J. Performance Data - General**

If IVD device, select "N/A" and the below criteria will be omitted from checklist. Performance data criteria relating to IVD devices will be addressed in Section K.

36) Full test report is provided for each completed test. A full test report includes: objective of the test, description of the test methods and procedures, study endpoint(s), pre-defined pass/fail criteria, results summary, conclusions, and an explanation of how the data generated from the test supports a finding of substantial equivalence.			X	
37)				
a) If there are requirements regarding performance data, such as special controls, in a device-specific regulation that are applicable to the device, the submission includes performance data to establish that the submitter has followed the device-specific requirement.			X	
b) If there is a device-specific guidance, other than a special controls guidance document, applicable to the device, the submission includes performance data to establish that the submitter has addressed the recommendations or otherwise has met the applicable statutory or regulatory criteria through an alternative approach.			X	
c) If there is a special controls document applicable to the device, the submission includes performance data to establish that the submitter has complied with the particular mitigation measures set forth in the special controls document or uses alternative mitigation measures but provides a rationale to demonstrate that those alternative measures identified by the firm will provide at least an equivalent assurance of safety and effectiveness.	X			
38) If literature is referenced in the submission, submission includes:			X	
39) For each completed nonclinical (i.e., animal) study conducted			X	

**K. Performance Characteristics - In Vitro Diagnostic Devices Only**

(Also see 21 CFR 809.10(b)(12))

Submission states that the device: (one of the below must be checked)

<input type="checkbox"/> is an in vitro diagnostic device.
<input checked="" type="checkbox"/> is not an in vitro diagnostic device.



**Decision:**  Accept  Refuse to Accept

If Accept, notify applicant.

If Refuse to Accept, notify applicant in writing and include a copy of this checklist.

**Digital Signature Concurrence Table**

Reviewer Sign-Off	Andrew I. Steen -S 2015.06.16 09:33:24 -04'00'
Branch Chief Sign-Off (digital signature optional)*	
Division Sign-Off (digital signature optional)*	

\* Branch and Division review of checklist and concurrence with decision required.  
Branch and Division digital signature optional.

## MEMO OF SOFTWARE REVIEW

**Name of the Device:** 3Shape Abutment Designer software  
**Common Name:** Abutment Designer  
**Applicant:** 3Shape A/S  
Copenhagen  
**Device Premarket Path:** K151455  
**Consult:** CON1512269  
**Lead Reviewer:** Andrew Steen  
**Kind of Device:** 21 CFR 872.3630; Endosseous dental implant abutment  
Class II. NNH  
**Predicate:** K100152 Sirona Dental CAD/CAM System  
**Date Sent:** 7/22/2015  
**Reviewer:** Catherine Li (CDRH, General Hospital Devices Branch),  
301-796-6304, [catherine.li@fda.hhs.gov](mailto:catherine.li@fda.hhs.gov)

(b)(4)



(b)(4)



## Software Architecture

(b)(4)



(b)(4)



**Original Submission Deficiency**


(b)(4)



(b)(4)



(b)(4)



Catherine Li  
CDRH/DAGRID/GHDB  
U.S. Food and Drug Administration  
Tel: (301) 796-6304

Digital Signature Concurrence Table	
Reviewer Sign-Off	<b>Catherine Li -A</b> <small>Digitally signed by Catherine Li-A DN: c=US, o=U.S. Government, ou=HHS, ou=FDA, ou=People, cn=Catherine Li-A, 0.9.2342.19200300.100.1.1=2000557465 Date: 2015.07.22 00:47:08 -04'00'</small>
Branch Chief Sign-Off	



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration  
 CDRH/ODE/DAGRID/DEDB  
 WO66 RM2604  
 10903 New Hampshire Ave  
 Silver Spring, MD 20993-0002  
 301-796-6284

**Premarket Notification [510(k)] Review**

<b>Date:</b> July 31, 2015			
<b>To:</b> FILE			
<b>From:</b> Andrew I. Steen, Mechanical Engineer			
<b>Subject:</b> Traditional 510(k)# K151455			
<b>Applicant:</b> 3Shape A/S		<b>Device Trade Name:</b> 3Shape Abutment Designer Software	
<b>Contact:</b> Hanne Nielsen		<b>Contact Title:</b> Regulatory Affairs Manager	
<b>Correspondent Firm:</b> 3shape A/S		<b>Phone:</b> 45 (702) 726-20 <b>Email:</b> hanne.nielsen@3shape.com	
<b>FDA Received Date:</b> June 1, 2015		<b>Due Date:</b> August 30, 2015	
<b>Reg #:</b> 872.3630 <b>Reg Name:</b> Endosseous Dental Implant Abutment		<b>Class:</b> II <b>Product Code(s):</b> NHA	
<b>Predicate Devices:</b>			
<b>Submission #</b>	<b>Pro Code</b>	<b>Device Trade Name</b>	<b>Owner</b>
K100152	NHA	Sirona Cad/cam System	Sirona Dental Systems GmbH
<b>Review Summary</b>			
<p>The subject device is a Endosseous Dental Implant Abutment with the following Indications for Use: "The 3Shape Abutment Designer Software is intended for designing endosseous dental implant abutments based on digital scans of 3D models. Intended users are Dental Practitioners and Dental Laboratory Staff. Intended Operational Environment is Dental Laboratories." It is for Rx use.</p>			
<b>Recommendation</b>			
<p>I recommend that the 3Shape Abutment Designer Software is/are in need of Additional Information (AINN)</p>			

**Review Team**

Lead Reviewer

Andrew I. Steen, Mechanical Engineer (CDRH/ODE/DAGRID/DEDB)



**I. Purpose and History**

3Shape A/S, of Copenhagen, Denmark, has submitted a Premarket Notification (510(k)) to introduce into U.S. interstate commerce *3Shape Abutment Designer™ Software*, a dental implant abutment type of device.

The original submission, received on June 1, 2015, 2013, consists of a cover letter, FDA submission cover sheets, an Indications for Use Statement (IFUS), a table of contents, a 510(k) Summary, a Truthful and Accurate Statement, a device description, a substantial equivalence comparison, a biocompatibility justification, EMC and electrical safety justification, shelf life justification, sterilization justification, corrosion testing justification, implant to abutment compatibility justification, material composition justification, deficiencies from previous submission K133457 and applicant's responses, applicable standards and FDA guidance, Standards Data Report forms (SDRs), software description, software verification and validation information, marketing brochure, and a device use manual.

The submission identifies prior submission K133457, CPT1200130, and I120810.

[TPLC Information](#) [Recall Information](#) [Historyfalls](#)

**II. Device/System Description**

Device Characteristics	Yes	No	Inadequate Or Marked
Is the intended use or fundamental technology new?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Is the device <u>life-supporting</u> or <u>life sustaining</u> ?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Are there any <u>direct</u> or <u>indirect</u> patient contacting components?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
• Is the device or a component an <u>implant</u> ?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Does the device use software/firmware?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
• Is the device, or does it contain, a <u>Mobile Medical App</u> ?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Does the device or a component need sterilization (by manufacturer or user)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
The device/system uses or is	a single use device(s) (SUD)		<input type="checkbox"/>
Is the device a <u>combination product</u> ?	N - Not a Part 3 Combination Product		<input type="checkbox"/>
Is the device electrical (battery or wall powered)?	Yes, it is mains powered Only		<input type="checkbox"/>
Check the attributes that are applicable to this submission.			
	Nanotechnology	Reprocessed SUD	Companion Diagnostic
Yes	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
No	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Unknown	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Device Description Table: Summary of important device characteristics			


(b)(4)



(b)(4)



(b)(4)



**III. Comparison of Indications for Use to Predicate Devices**

Comparison of Indications for Use								
<b>Subject</b>								
510(k) #: K151455								
	Rx/OTC: Rx							
Intended Population	Adults Only	Adults and Pediatrics	Transitional Adolescent A	Transitional Adolescent B	Adolescent	Child	Infant	Neonate/ Newborn
Yes	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
No	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Unknown	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

**Comparison of Indications for Use**

Indications for Use: "The 3Shape Abutment Designer Software is intended for designing endosseous dental implant abutments based on digital scans of 3D models. Intended users are Dental Practitioners and Dental Laboratory Staff. Intended Operational Environment is Dental Laboratories."

**Predicate(s)**

510(k)#: K100152

Rx/OTC: Rx

Intended Population:

Indications for Use: The sirona dental cad/cam system is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations. The system consists of three major parts: tibase, incoris mesostructure, and cad/cam software. Specifically, the incoris mesostructure and tibase components make up a two-piece abutment which is used in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity. The incoris mesostructure may also be used in conjunction with the camlog titanium base cad/cam (types k2244.Xxxx) k083496) in the camlog implant system. The cad/cam software is intended to design and fabricate the incoris mesostructure. The incoris mesostructure and tibase two-piece abutment is compatible with the following implants system:

- nobel biocare replace (k020646)
- nobel biocare branemark (k022562)
- fridament xive (k013876)
- biomet 3i osseotite (k980549)
- astra tech osseospeed (k091239)
- zimmer tapered screw-vent (k061410)
- straumann synocta (k061176)

Indications for Use Table: Compares the indications for use of the subject and predicate devices.

**Reviewer Recommendation**

The Comparison of the Indications for Use is [not] acceptable.

**IV. Comparison of Technology to Predicate Devices**

<b>Device &amp; Predicate Device(s):</b>	<u>K151455</u>	<u>K100152</u>
<b>General Device Characteristics</b>		
Trade Name		
Applicant Company		

Predicate is a scanner, milling unit, design software, and a two-piece abutment.

Submission device is CAD/CAM only

**Reviewer Recommendation**

The Comparison of the Technology to Predicate Devices is [not] acceptable.

**V. Labeling**

Labeling Review Needed?	Yes <input checked="" type="checkbox"/>	No <input type="checkbox"/>
Usability Consult Needed?	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>

**A General Labeling Requirements**

General Labeling Requirements	Yes	No N/A	Inadequate or Marked
Is the prescription statement (or "Rx only") included?	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
The indications for use are consistent with the IFU page?	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Appropriate contraindications provided?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Appropriate warnings provided?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Instructions are in accordance with the guidance (if applicable)?	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Appropriate labeling inside device?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Appropriate label/indicator outside device?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Appropriate Manual labeling?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
What MRI safety information does the labeling contain?	Not Needed		<input type="checkbox"/>
Labeling Table: A summary of the adequacy of several labeling requirements.			

**Reviewer Recommendation**  
 (b)(4)

**VI. Cleaning, Disinfection, Sterilization, Shelf-Life and Reuse**

Sterility Review Needed?	Yes	<input type="checkbox"/>	No	<input checked="" type="checkbox"/>
Sterility Consult Needed?	Yes	<input type="checkbox"/>	No	<input checked="" type="checkbox"/>

The submission states that the device is only software and therefore shelf life and sterilization requirements do not apply.

**Reviewer Recommendation**  
 Cleaning, Sterilization, Shelf-Life and Reuse descriptions are acceptable.

**VII. Biocompatibility**

Biocompatibility Review Needed?	Yes	<input type="checkbox"/>	No	<input checked="" type="checkbox"/>
Biocompatibility Consult Needed?	Yes	<input type="checkbox"/>	No	<input checked="" type="checkbox"/>


The submission states that the device has no patient contacting components and thus biocompatibility requirements are not apply. The submission states that identification of the material is also not applicable.

**Reviewer Recommendation**  
 The Biocompatibility is acceptable.

**VIII. Software/Firmware**

Software Review Needed?	Yes	<input checked="" type="checkbox"/>	No	<input type="checkbox"/>
Software Consult Needed?	Yes	<input type="checkbox"/>	No	<input checked="" type="checkbox"/>

(b)(4)



**IX. EMC & Electrical, Mechanical and Thermal Safety & Risk Analysis**

EMC Review Needed?	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>
EMC Consult Needed?	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>


The submission states that the submission device has not patient or user contacting components as it is just software. The submission states that EMC and Electrical safety evaluation is not applicable.

**Reviewer Recommendation**

The EMC, EMT and Risk Analysis are acceptable.

**X. Performance Testing**

(b)(4)



**XI. Kit Certification**

This section is N/A.

**XII. Manufacturing Information**

This section is N/A.

**XIII. References**

The submission references the following standards or guidance documents and has provided the appropriate Standards Data Report Forms.


- Class II Special Controls Guidance Document: Root-form Endosseous Dental Implants and Endosseous Dental Implants
- ISO 13485:2015 Quality management systems – requirements for regulatory process
- ISO 14971:2012 Application of Risk Management to Medical Devices
- IEC 62304:2012 Software life cycle processes

**XIV. SE Flowchart Questions**

Substantial Equivalence Determination	Yes	No
Is the predicate device legally marketed?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Do the devices have the same intended use?	<input type="checkbox"/>	<input type="checkbox"/>

**XV. Original Deficiencies**

(b)(4)





(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



<b>Digital Signature Concurrence Table</b>	
Reviewer Sign-Off	Andrew I. Steen - S 2015.07.31 12:25:23 -04'00'

XVI. Contact History



## MEMO OF SOFTWARE REVIEW

**Name of the Device:** 3Shape Abutment Designer software  
**Common Name:** Abutment Designer  
**Applicant:** 3Shape A/S  
Copenhagen  
**Device Premarket Path:** K151455/S001 Traditional  
**Consult:** CON162977  
**Lead Reviewer:** Andrew Steen  
**Kind of Device:** 21 CFR 872.3630; Endosseous dental implant abutment  
Class II. NNH  
**Predicate:** K100152 Sirona Dental CAD/CAM System  
**Date Sent:** 1/19/2016  
**Reviewer:** Catherine Li (CDRH, General Hospital Devices Branch),  
301-796-6304, [catherine.li@fda.hhs.gov](mailto:catherine.li@fda.hhs.gov)

(b)(4)



(b)(4)



(b)(4)



†

(b)(4)



(b)(4)



(b)(4)



**Original Submission Deficiency**

(b)(4)



(b)(4)



(b)(4)





**Deficiencies:**

(b)(4)



Catherine Li

CDRH/DAGRID  
U.S. Food and Drug Administration  
Tel: (301) 796-6304

Digital Signature Concurrence Table	
Reviewer Sign-Off	Catherine Li -A <small>Digitally signed by Catherine Li -A DN: cn=US, o=U.S. Government, ou=HHS, ou=FDA, ou=People, c=Catherine Li -A, e=Catherine.Li@FDA.HHS.gov, o=U.S. Government, ou=FDA, ou=People, ou=Catherine Li -A, o=Catherine Li -A, o=Catherine Li -A, o=Catherine Li -A Date: 2016.02.19 12:23:26 -0500</small>
Branch Chief Sign-Off	

## MEMO OF SOFTWARE REVIEW

**Name of the Device:** 3Shape Abutment Designer software  
**Common Name:** Abutment Designer  
**Applicant:** 3Shape A/S  
Copenhagen  
**Device Premarket Path:** K151455/S001 Traditional  
**Consult:** CON162977  
**Lead Reviewer:** Andrew Steen  
**Kind of Device:** 21 CFR 872.3630; Endosseous dental implant abutment  
Class II. NNH  
**Predicate:** K100152 Sirona Dental CAD/CAM System  
**Date Sent:** 1/19/2016, updated on 3/8/2016  
**Reviewer:** Catherine Li (CDRH, General Hospital Devices Branch),  
301-796-6304, [catherine.li@fda.hhs.gov](mailto:catherine.li@fda.hhs.gov)

(b)(4)



(b)(4)



(b)(4)



†

(b)(4)



(b)(4)



(b)(4)

A large black rectangular redaction box covering the majority of the upper half of the page.

Original Submission Deficiency

(b)(4)

A very large black rectangular redaction box covering the entire middle and lower portion of the page.



(b)(4)



(b)(4)



**Deficiencies:**


(b)(4)



(b)(4)



(b)(4)



Catherine Li  
CDRH/DAGRID  
U.S. Food and Drug Administration  
Tel: (301) 796-6304

Digital Signature Concurrence Table	
Reviewer Sign-Off	Catherine Li -A <small>Digitally signed by Catherine Li -A DN: c=US, o=U.S. Government, ou=HHS, ou=FDA, ou=People, cn=Catherine Li -A, 09.2342.19208300.100.1.1=2000657465 Date: 2016.03.08 11:57:06 -05'00'</small>
Branch Chief Sign-Off	



Food and Drug Administration  
 CDRH/ODE/DAGRID/DEDB  
 WO66 RM2604  
 10903 New Hampshire Ave  
 Silver Spring, MD 20993-0002  
 301-796-6284

### Premarket Notification [510(k)] Review

**Date:** August 18, 2016  
**To:** FILE  
**From:** Andrew I. Steen, Mechanical Engineer  
**Subject:** Traditional 510(k)# K151455/S001

<b>Applicant:</b> 3Shape A/S	<b>Device Trade Name:</b> 3Shape Abutment Designer Software
<b>Contact:</b> Hanne Nielsen	<b>Contact Title:</b> Regulatory Affairs Manager
<b>Correspondent Firm:</b> 3shape A/S	<b>Phone:</b> 45 (702) 726-20 <b>Email:</b> hanne.nielsen@3shape.com
<b>FDA Received Date:</b> January 11, 2016	<b>Due Date:</b> February 10, 2016
<b>Reg #:</b> 872.3630 <b>Reg Name:</b> Endosseous Dental Implant Abutment	<b>Class:</b> II <b>Product Code(s):</b> PNP

#### **Predicate Devices:**

Submission #	Pro Code	Device Trade Name	Owner
K100152	NHA	Sirona Cad/cam System	Sirona Dental Systems Gmbh

#### **Review Summary**

The subject device is an Endosseous Dental Implant Abutment with the following Indications for Use: "The 3Shape Abutment Designer Software is intended for designing endosseous dental implant abutments based on digital scans of 3D models. Intended users are Dental Practitioners and Dental Laboratory Staff. Intended Operational Environment is Dental Laboratories." It is for Rx use.

#### **Recommendation**

I recommend that the 3Shape Abutment Designer Software is/are **Substantially Equivalent (SESE)**

#### Review Team

Lead Reviewer  
 Software Consultant

Andrew I. Steen, Mechanical Engineer (CDRH/ODE/DAGRID/DEDB)  
 Catherine Li, Electrical Engineer (CDRH/ODE/DAGRID/GHDB)

(b)(4)



(b)(4)





(b)(4)



(b)(4)



(b)(4)



(b)(4)




(b)(4)



(b)(4)

II



(b)(4)



(b)(4)



v



(b)(4)



(b)(4)



(b)(4)

**VII. Biocompatibility**

Biocompatibility Review Needed?	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>
Biocompatibility Consult Needed?	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>

The submission states that the device has no patient contacting components and thus biocompatibility requirements are not apply. The submission states that identification of the material is also not applicable.

**Reviewer Recommendation**

The Biocompatibility is acceptable.

**VIII. Software/Firmware**

Software Review Needed?	Yes <input checked="" type="checkbox"/>	No <input type="checkbox"/>
Software Consult Needed?	Yes <input checked="" type="checkbox"/>	No <input type="checkbox"/>

(b)(4)

(b)(4)



**IX. EMC & Electrical, Mechanical and Thermal Safety & Risk Analysis**

EMC Review Needed?	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>
EMC Consult Needed?	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>

(b)(4)

**Reviewer Recommendation**

The EMC, EMT and Risk Analysis are acceptable.

**X. Performance Testing**

(b)(4)

**Reviewer Recommendation**

The Performance Testing [Verification & Validation] is acceptable.

**XI. Kit Certification**

This section is N/A.

**XII. Manufacturing Information**

This section is N/A.

**XIII. References**

The submission references the following standards or guidance documents and has provided the appropriate Standards Data Report Forms.

- ISO 13485:2015 Quality management systems – requirements for regulatory process
- ISO 14971:2012 Application of Risk Management to Medical Devices
- IEC 62304:2012 Software life cycle processes

**XIV. SE Flowchart Questions**

Substantial Equivalence Determination

Yes

No

(b)(4)



**XV. Original Deficiencies**

(b)(4)



(b)(4)



(b)(4)





(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)





(b)(4)



**XVI. April 23 Interactive Deficiencies**

(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)





(b)(4)



(b)(4)



**XVII. June 22 Interactive Deficiencies**

(b)(4)



(b)(4)



(b)(4)



(b)(4)



**XVIII. Contact History**

(b)(4)



(b)(4)



Digital Signature Concurrence Table	
Reviewer Sign-Off	Andrew I. Steen -S 2016.08.18 16:01:08 -04'00'

**Online Implant Libraries Provider List as of May 25, 2016**







3Shape Headquarter Europe, Middle East & Africa  
 Holmens Kanal 7  
 1060 Copenhagen, Denmark  
 Tel: +45 7027 2620

3Shape Asia  
 Room 906, Tower A of Eton Place  
 No. 69, Dongfang Road  
 200120 Shanghai, China  
 Tel: +86 21 5835 2281

3Shape Latin America  
 Carrera 13 # 82-91  
 Oficina 401  
 110221 Bogotá, Colombia  
 Tel: +57 1691 9508

3Shape North America  
 Somerset Hills Corporate Center  
 10 Independence Boulevard, Suite 150  
 Warren, New Jersey 07059, USA  
 Tel: +1 (908) 867 0144

More information at  
 www.3shape.com  
 info@3shape.com

Implant systems	Providers of implant solutions	Implant systems	Implant Solutions	Manu- facturing	Consumables (for local use)	Additional info
Heraeus Kulzer (DE) www.heraeus.com	Thornion	X	X	X		
Implant Direct (US)* www.implantdirect.com	Biorhance	X	X	X		FDA CE, Health Canada
Intra-Lock International (US)* www.intra-lock.com	INTRA-LOCK	X	X	X		FDA CE, Health Canada
Vocalar Wierand (DE) www.wierand-dental.de	Constaris Implants	X	X	X		
Lasak (CZ) www.lasak.com	Lasak Bridge	X	X	X		CE, ISO 9001
LeStructura (IT) www.lestructura.it	Bridges & Maxilla - Implant - and many more	X	X	X		CE, ISO 9001
Meditentika (DE) www.meditentika.de	Meditentika M- Implant	X	X	X		CE
Meditentika Medical (DE) www.meditentika.de	ECT Implant	X	X	X		CE
Medical Instinct (DE) www.medical-instinct.de	Medical Instinct	X	X	X		CE
Medical Production (FR) www.medical-production.eu	EuroMax	X	X	X		CE, ISO 9001
Megagen Implant (KR)* www.megagen.com	AnyAge Maxilla - AnyAge Maxilla - DCS Level	X	X	X		FDA, CE, ISO 9001
MIS Implants Technology (IL)* www.mis-implants.com	MIS Implants	X	X	X		CE, FDA
Nedent (BR) www.nedent.com.br	Nedent Implant Systems	X	X	X		
Necos (UK)* www.necos.com	Necos	X	X	X		CE, FDA
Nobel Biocare (CH)* www.nobelbiocare.com		X	X	X		CE, FDA
NIT Trading (DE)* www.nit-trading.com	Thornion - Omega & Maxilla - IT	X	X	X		CE, FDA, CAD/CAM, ISO 9001, TQC, Mexico Certification
Phibo (ES) www.phibo.com	Phibo - and many more	X	X	X		CE, ISO 9001
Provalat (DE) www.provalat.de	Provalat	X	X	X		CE, EN ISO 9001:2008, EN ISO 13485:2003, CE, CE Marking Directive 93/42/EWG
Ritter Implants (DE) www.ritterimplants.com	Age-Plus Tec - MIS Implants	X	X	X		CE
Straumann (CH)* www.straumann.com		X	X	X		CE, FDA



More information at  
 www.3shape.com  
 info@3shape.com

3Shape North America  
 Somerset Hills Corporate Center  
 10 Independence Boulevard, Suite 150  
 Warren, New Jersey 07059, USA  
 Tel: +1 (908) 867 0144

3Shape Latin America  
 Carrera 13 # 82-91  
 Oficina 401  
 110221 Bogotá, Colombia  
 Tel: +57 1691 9508

3Shape Asia  
 Room 906, Tower A of Eton Place  
 200120 Shanghai, China  
 Tel: +86 21 5835 2281

3Shape Headquarters  
 Europe, Middle East & Africa  
 Holmens Kanal 7  
 1060 Copenhagen, Denmark  
 Tel: +45 7027 2620

- Only Implant Libraries marked with an asterisk (\*) are approved for use in the United States. See the full list of approved Implant Systems in the Appendix.
- The above information is solely provided by third parties, and 3Shape assumes no responsibility for its accuracy or completeness.

Please note:

Implant systems	Providers of implant solutions	Additional Info
Sweden&Martina S.p.A (IT) www.sweden-martina.com	Sweden&Martina	
Target3D (USA) www.target3d.com	Osstem - Biocrown - Biocrown - and more	- CE
Thommen Medical (CH)* www.thomme-medical.com	Thommen	- CE, FDA
TRI Implants (CH) www.tri-implants.com	TRI Implants	
Vulcan Custom Dental (US)* www.vulcandental.com	Biocrown	- FDA
Zfx (DE) www.zfx-dental.com	Biogyrus - Titanium - Ceramic - and more	- CE
Zimmer Dental (US)* www.zimmerdental.com		- FDA - CE

Implant systems	Providers of implant solutions	Additional Info
Nobel Biocare*		
Bionix 3i*		
Streamline*		
Zimmer Dental*		
DENTSPLY - Aera Tech*		
DENTSPLY - Ribdent*		
CAMLOG*		
Other systems		
1- Piece (Titanium)		
1- Piece (Zirconia)		
2- Piece (Ti base + Zr abutment)		
Implant Bars & Bridges		
Digital Implant models		
Providers milling center		
Local milling by 3rd party		
Bands with pre-milled interface (for 2-piece abutment)		
Titanium bases (for 2-piece abutment)		
Interoral Scan bodies		
Digital model analogs		
Multi-unit abutments		
Consumables (for local use)		
Manu-facturing		
Implant Solutions		
Comments		



More information at  
www.sshape.com  
info@sshape.com

3Shape North America  
Somerset Hills Corporate Center  
10 Independence Boulevard, Suite 150  
Warren, New Jersey 07059, USA  
Tel: +1 (908) 867 0144

3Shape Latin America  
Carretera 13 # 82-91  
Oficina 401  
110221 Bogotá, Colombia  
Tel: +57 1691 9508

3Shape Asia  
Room 906, Tower A of Eton Place  
No. 69, Dongfang Road  
200120 Shanghai, China  
Tel: +86 21 5835 2281

3Shape Headquarters  
Europe, Middle East & Africa  
Holmens Kanal 7  
1060 Copenhagen, Denmark  
Tel: +45 7027 2620

Implant System	Implant System	Implant System	Implant System
AB-Dental A.B. Dental Implant System K13125	Alpha - Bio Tec Alpha-Bio Tec K06364	Anthogyr AXIOM <sup>®</sup> REG K131066 AXIOM <sup>®</sup> Dental Implant System K101913 AXIOM <sup>®</sup> 2.8 K141450	DYNA Dental DYNA Implant Design K930868
AVINENT AVINENT Dental Implant System K121873	ATLANTIS ATLANTIS ABUTMENT AND ATLANTIS ABUTMENT SCREW K981858	BTI BTI Dental Implant 5.5-6.5 K091387 BTI Dental Implant Tiny <sup>®</sup> 3.0 K092112 BTI Internal Dental Implant System K053555	BEGO Implant Systems BEGO Semados <sup>®</sup> S-Line K090716
Bicon Bicon Dental Implant System 3.0mm K101849 Bicon Implants with a 2.5mm Internal Connection K092035 Bicon 5.0 x 5.0mm and 6.0 x 5.0 dental Implants K073368 Bicon 4.5, 6.0mm and 6.0*6.0 mm dental Implants K062044 Bicon 5.0 x 6.0mm dental Implant K050408 Bicon Surface Treatment K042637 Bicon the 5.0 x 6.0mm dental implant K023705 Bicon Transitional Implant (TRI) system K010185 Bicon 6.0 x 5.7mm dental implant K994037 Bicon the 4.5mm diameter bicon dental implant K982488 Bicon bone screw system K972417 Bicon Dental Implants II part abutment system K972029	Biodenta Biodenta Dental Implant System - Bone Level Tapered D.3.0 K133884 Biodenta Dental Implant System - Bone Level D.3.0 to 6.0 mm K123512 Biodenta Dental Implant System - Bone Level Tapered K123415 Biodenta Dental Implant System - Bone Level Tapered K111003 Biodenta Dental Implant System K093630	BioHorizons BioHorizons Tapered Internal Implants K143022 BioHorizons Tapered Internal Plus Implants K121787 BioHorizons Laser-Lok 3.0 Implant System K093321 BioLok Micro-Lok Implant System K081818 BioHorizons Internal Implant System K073282 BioHorizons Internal Implant System K073268 BioHorizons Tapered Internal Implant System K071638 BioHorizons Single-Stage Implant K053152 Maximus <sup>™</sup> 3.0mm and 05 Implants K052419 Prodigy System <sup>™</sup> Dental Implants K042429	Biomet 3i 3iTX <sup>®</sup> External Hex Dental Implants K133049 CP4 OSSEOTITE <sup>™</sup> Certain <sup>™</sup> Dental Implants K130549 3iTX Dental Implant K122300 OSSEOTITE <sup>™</sup> 2 Dental Implants K111216 Full OSSEOTITE <sup>™</sup> Certain <sup>™</sup> II Dental Implant K100724 BIOMET 3i Nanotite <sup>™</sup> Dental Implants K072363 Certain <sup>™</sup> PREVAL <sup>™</sup> Implants K061629 OSSEOTITE NT <sup>™</sup> Implant System K014235 Biomet 3i OSSEOTITE <sup>™</sup> K980549
CAMLOG CAMLOG SCREW-LINE K083496 CAMLOG ROOT-LINE (I-Line) K000100 CONELOG SCREW-LINE K13779 Sy K133991 CAMLOG CYLINDER K000065 CAMLOG SCREW K000099	Cortex Dental Implants SATURN Dental Implant System K131258 CORTEX Dental Implant System K090709	Credent Credent Solides Customized Abutment K150012 Credent Solides Custom Abutment K113738	Dentium Dentium Slim Onebody System K11162 Dentium S Implants II Model S0k483408r, S0k483 K102308 Dentium Speedy Orthodontic Screw K060500 Dentium Implantum II K060501 Dentium Co., Ltd. Implantum K041368
DENTISPLY - Astra Tech Osseospeed <sup>™</sup> Profile EV K130999 Osseospeed <sup>™</sup> Plus K120414 Astra Tech Implant System Plus K111287 Astra Tech Implant System K101732 Astra Tech Implant System K091239 Osseospeed <sup>™</sup> Profile System K080156 Osseospeed <sup>™</sup> Narrow K080396 Osseospeed <sup>™</sup> 4.05-6mm K063779 Astra Tech Implants-Dental System K041492	DENTISPLY - Friadent ANKRYOS <sup>®</sup> C/X Implant System K140347 ANKRYOS <sup>®</sup> C/X Dental Implant System K083805 FRIALIT <sup>®</sup> plus Dental Implant System, XIVE <sup>®</sup> 5 plus Dental Implant System, K073075 ANKRYOS <sup>®</sup> Dental Implant System / Claims K073067 ANKRYOS <sup>®</sup> Dental Implant System K041509 FRIADENT plus Dental Implant Systems K040170 ANKRYOS <sup>®</sup> Dental Implant System K040946 XIVE <sup>®</sup> Dental Implant System K013867 ANKRYOS <sup>®</sup> Dental Implant System K012681		

The following implant systems have been approved for use in the USA by the FDA. Notice, the approvals apply only for the implant systems and not for the implant libraries above. For further information please refer to the 510(k) numbers listed below with the implant systems.

Appendix: Regulatory Declaration



Product Name	Device ID	Manufacturer
DIO DIO UP HSA Internal Sub-Merged Implant System	K122519	DIO
DIO DIO STEADY External Implant System	K112746	DIO
DIO DIO STEADY External Implant System (2.5/3.0 mm)	K100100	DIO
SM-Extra Wide (BBM) Implant System	K080128	DIO
DIO Protein Implant System	K080126	DIO
DIO BioTite-H Implant System	K073070	DIO
DIO Implant System	K070570	DIO
SM <sup>®</sup> Internal/External Implant System	K070569	DIO
DIO Protein Implant System	K070568	DIO
SM <sup>®</sup> Implant Systems	K061797	DIO
<b>EurotekniKa</b>		
UNIVERSAL, AESTHETICA, MATEA, NATURA	K083670	
Obi	K083669	
<b>GC America</b>		
GC Aadvio Implant System	K093749	
JNE Implant System	K072425	
<b>Green DenTech</b>		
DENRACLE DENTAL ABUTMENT	K131468	
<b>ImplantDirect</b>		
Spectra System	K061319	
IntroActive/SwissPlus2 Implant System	K130572	
Spectra System Dental Implants 2008	K090234	
Swissplant Dental Implant System	K081396	
<b>Intra-Lock</b>		
Intra-Lock <sup>®</sup> OP Dental Implants	K130140	
Intra-Lock <sup>®</sup> Dental Implants	K133613	
Intra-Lock Dental Implant System with Blossom	K103194	
Mini Drive-Lock <sup>™</sup> Dental Implant System	K070601	
Intra-Lock MIILO <sup>™</sup> Dental Implant System	K050970	
Intra-Lock Hydroxyapatite Coated Implants	K031322	
Intra-Lock Transitional Implant System	K021915	
Intra-Lock Dental Implant System	K021322	
<b>Keystone Dental</b>		
Genesis Implant System	K101545	
Olympus Dental Implant System	K071070	
<b>Lifecore Biomedical</b>		
Lifecore FirmoConnect <sup>™</sup> Internal Connection Implant System	K051614	
Lifecore FirmoSolo <sup>™</sup> One-Piece Implant System	K050506	
RENOVA <sup>™</sup> Internal Hex Implant System	K032774	
Restore <sup>™</sup> RBM Self-Tapping Regular Diameter Dental Implant	K002037	
Lifecore Stage-1 Single Stage RBM Dental Implant System	K002226	
Single Stage Dental Implant	K091114	
Stage-1 RDS TFS Dental Implant System	K094205	
<b>MIS Implants</b>		
Conical Connection Implants	K112162	
SEVEN Implants, BioCom Implants, LANCE Implants	K103089	
UNO Narrow Implant	K092555	
UNO -One Piece Screw-Type Dental Implant	K080162	
Misral One Stage Dental Implant	K070022	
MIS Dental Implant System	K040807	
<b>Megagen</b>		
AnyOne <sup>™</sup> Internal Implant System	K123988	
Xped AnyBridge Internal Implant System	K122231	
Xped AnyBridge Internal Implant System	K123870	
RESQ <sup>™</sup> External Implant System	K081302	
RESQ <sup>™</sup> Internal Dental Implant System	K073058	
RESQ <sup>™</sup> Internal Dental Implant System	K062166	
RESQ <sup>™</sup> Implant System	K053333	
Intermezzo <sup>™</sup> Plus	K053354	
ExFect <sup>™</sup> Implant Systems	K052369	
Intermezzo <sup>™</sup> Implant Systems	K051018	
<b>Ness</b>		
Ness Proactive Tapered Implant	K113376	
Ness Implant System 3.25	K090452	
Ness Proactive Implant	K083561	
Ness Implant System	K041395	
Ness Coccr Abutments	K150669	
Ness Ti Reinforced Membrane	K143327	
Ness Access Abutments	K081851	
Ness Various Titanium Abutments	K071838	
<b>DIO</b>		
Nobel Biocare		
NobelActive Wide Platform (WP)	K133731	
NobelActive 3.0	K102436	
NobelActive 8.5mm & 18.0mm	K083205	
NobelReplace Hexagonal Implant	K073142	
NobelPerfect Conical Connection Hexagonal	K072980	
NobelActive Internal Connection Implant	K071370	
NobelReplace Tapered Conical Connection	K062566	
NobelReplace	K022562	
NobelReplace	K020646	
<b>Ostem</b>		
M5 SA Implant System	K122171	
T5 Fixture System	K121995	
T5 Implant System	K121585	
ET/55 Implant System	K120847	
HG II Short Fixture System	K091678	
GS III System	K091208	
M5 System	K083067	
H5 II Short Fixture System	K083633	
GS III Fixture System	K082213	
NT-IT System	K081078	
HG II Fixture System	K080744	
M5 System (Narrow Ridge)	K080594	
HU II - HS II Fixture System	K080387	
US-G5 System	K073247	
US-G5 Ultra Wide System	K073465	
M5 System (Denture)	K072959	
GS Fixture System	K072896	
GS System	K063861	
US System	K062030	
55 System	K062051	
AVANA Dental Implant System	K051576	
<b>PrismaK DentalCraf</b>		
Inclusive <sup>™</sup> Tapered Implant System	K121406	
Inclusive <sup>™</sup> Mini Implant	K100932	
<b>Ritter Implants</b>		
Ritter Implant System	K131557	
<b>Straumann</b>		
Straumann SLActive <sup>®</sup> and Roxolid <sup>®</sup>	K130222	
SLA <sup>®</sup> , SLActive <sup>®</sup> and Roxolid <sup>®</sup> Product Families	K12784	
Straumann Bone Level 04.1 mm and 04.8 mm Roxolid	K122855	
Straumann Bone Level (BL) 04.1 mm and 04.8 mm Regular Connection (RC) Roxolid	K121313	
Straumann Narrow Neck CrossFit <sup>®</sup> (NNC) 03.3 mm Dental Implant System	K11357	
Straumann Dental Implant System	K083550	
Straumann Dental Implant System	K081419	
Straumann Bone Level	K062129	
Straumann synOcta <sup>®</sup>	K061176	
Straumann Narrow Neck Implants	K060958	
Straumann SLActive <sup>®</sup> Implants	K033088	
ITI Dental Implant System	K033984	
Straumann <sup>®</sup> Bone Level	K030007	
Straumann <sup>®</sup> Bone Level	K010291	
<b>SYBRON</b>		
Pitt-Easy Dental Implant System	K083297	
<b>Thommen Medical</b>		
SPI <sup>®</sup> Dental Implant, ELEMENT	K093615	
SPI <sup>®</sup> System Dental Implants	K090154	
SPI <sup>®</sup> CONTACT Platform 4.0 mm	K072933	
SPI <sup>®</sup> ELEMENT Platform 4.0 mm	K070007	
SPI <sup>®</sup> Dental Implant	K051502	
SPI <sup>®</sup> CONTACT Dental Implant	K034014	
<b>3Shape Latin America</b>		
Correa 13 # 82-91		
Oficina 401		
110221 Bogotá, Colombia		
Tel: +57 1691 9508		
<b>3Shape Asia</b>		
Room 906, Tower A of Eton Place		
No. 69, Dongfang Road		
200120 Shanghai, China		
Tel: +86 21 5835 2281		
<b>3Shape Europe, Middle East &amp; Africa</b>		
Holmens Kanal 7		
1060 Copenhagen, Denmark		
Tel: +45 7027 2620		
<b>3Shape North America</b>		
Somersset Hills Corporate Center		
10 Independence Boulevard, Suite 150		
Warren, New Jersey 07059, USA		
Tel: +1 (908) 867 0144		
www.3shape.com		
info@3shape.com		



More information at  
 www.3shape.com  
 info@3shape.com

3Shape North America  
 Somerset Hills Corporate Center  
 10 Independence Boulevard, Suite 150  
 Warren, New Jersey 07059, USA  
 Tel: +1 (908) 867 0144

3Shape Latin America  
 Carrera 13 # 82-91  
 Oficina 401  
 110221 Bogotá, Colombia  
 Tel: +57 1691 9508

3Shape Asia  
 Room 906, Tower A of Eton Place  
 No. 69, Dongfang Road  
 200120 Shanghai, China  
 Tel: +86 21 5835 2281

3Shape Headquarters  
 Europe, Middle East & Africa  
 Holmens Kanal 7  
 1060 Copenhagen, Denmark  
 Tel: +45 7027 2620

K133339	Zimmer Dental Tapered Screw-Vent® T Implant, HA Coated; Zimmer Dental Tapered
K132258	Zimmer Dental Tapered Screw-Vent® X Implant System
K113753	Zimmer Dental Tapered Screw-Vent® X Implant
K112160	Zimmer Dental Tapered Screw-Vent® X Implant
K111889	Zimmer Dental Tapered Screw-Vent® M Implant
K101880	Zimmer Dental Tapered Screw-Vent® T Implant;
K101977	Zimmer Dental Tapered Screw-Vent® T Implant;
K093164	3.25mm Spline® Twist™ Implant, HA Coated
K082639	Zimmer Dental Tapered Screw-Vent® Implants
K072589	Zimmer Dental Tapered Screw-Vent® Implant 4.1mm
K071235	Zimmer® One-Piece Implant, 3.0mm, Angled
K063523	Zimmer® One-Piece Implant, 3.7mm and 4.7mm
K061410	Zimmer Dental Tapered Screw-Vent® Implant System SV; Zimmer Dental Screw-Vent®
K061717	Zimmer® One-Piece Implant 3.7mm Angled
K062281	Zimmer® One-Piece Implant 4.7mm Straight
K052997	Zimmer® One-Piece Implant
K142082	Zimmer 3.1mm Dental Implant 2.9mm Angled



**List of Implant Providers and Implant Library information as included on the 3Shape US Server**

# Implant Libraries for the 3Shape Dental System™

June, 2016

Implant systems	Providers of implant solutions	Additional Info	Comments
Argen (US)	3Shape www.agen.com	X	
Atlantis™ (SE)	3Shape www.dentsplyimplants.com	Information not available	
Bicon (US)	3Shape www.bicon.com	X	
Biodenta (CH)	3Shape Biodenta Bone Level (no tissue level) Biodenta Custom (Custom) TSL Biodenta Bone Level and Lock - Boost	X	Implant bridges only available for Abutment level, not implant level. Biodenta Contact Manufacturer and comply to QSR
Biomet (US)	3Shape www.biomet.com	Information not available	
BlueskyBio (US)	BlueskyBio	X	
CAMLOG (CH)	3Shape www.camlog.com	X	
CAP (US)	3Shape www.cap-us.com	X	
CNC (US)	3Shape www.custom-milling.com	X	
Core3D International (NL)	3Shape www.core3denters.com	Information not available	
Cortex Dental (IL)	3Shape www.cortex-dental.com	X	
Credent (US)	3Shape www.credent.com	X	
Dentium (KR)	3Shape www.dentium.com	X	
Euroteknika (FR)	3Shape www.euroteknika.com	X	
Gildwell (US)	3Shape www.gildwell-dental.com	X	
Hiossen (US)	3Shape www.hiossen.com	Information not available	
Implant Direct (US)	3Shape www.implantdirect.com/cusom-direct	X	
Intra-Lock International (US)	3Shape www.intra-lock.com	X	
Medentika (DE)	3Shape www.medentika.de	X	
Megagen Implant (KR)	3Shape www.megagen.com	X	
3Shape Asia No. 69, Dongfang Road Room 906, Tower A of Eton Place Carrera 13 # 82-91 Somerset Hills Corporate Center 3Shape North America www.3shape.com More information at www.3shapedental.com	3Shape www.3shape.com	3Shape Asia No. 69, Dongfang Road Room 906, Tower A of Eton Place Carrera 13 # 82-91 Somerset Hills Corporate Center 3Shape North America www.3shape.com More information at www.3shapedental.com	3Shape Asia No. 69, Dongfang Road Room 906, Tower A of Eton Place Carrera 13 # 82-91 Somerset Hills Corporate Center 3Shape North America www.3shape.com More information at www.3shapedental.com



Implant systems	Implant systems							Implant Solutions					Manu- facturing		Consumables (for local use)					Additional Info		
	Nobel Biocare	Biomet 3i	Straumann	Zimmer Dental	DENTSPLY - Astra Tech	DENTSPLY - Frident	CAMLOG	Other systems	1-Piece (Titanium)	1-Piece (Zirconia)	2-Piece (Ti base + Zr abutment)	Implant Bars & Bridges	Digital Implant models	Providers' milling center	Local milling by 3rd party	Blanks with pre-milled interface (for 1-piece abutment)	Therium bases (for 2-piece abutments)	Intraoral Scan bodies	Digital model analogs		Multi-unit abutments	Cleared for use in US
MIS Implants Technologies Ltd. (IL) - FDA cleared only www.mis-implants.com							MIS Implants			X	X			X	X	X	X	X		X	X	-CE
Neodent (BR) www.neodent.com.br	👍		👍				Neodent Implant Systems	X	X	X	X				X			X			X	
Neoss (UK) - all libraries www.neoss.com	👍	👍	👍	👍	👍		Neoss	X		X	X	X	X	X	X	X	X	X	X	X	X	-CE
Nobel Biocare (CH) www.nobelbiocare.com	👍	👍	👍	👍	👍			X	X	X			X								X	-CE
NTTrading (DE) www.nt-trading.com	👍	👍	👍	👍	👍	👍	Thommen Sweden & Martina BEGO	X		X	X			X	X	X	X	X	X		X	-CE, CMDCAS, GOST R, TGC - Mexico Certification
Paltop (IL) paltopdental.com	Information not available																				X	
Straumann (CH) - FDA cleared only www.straumann.com			👍					X	X	X				X							X	-CE
Sweden&Martina S.p.A (IT) www.sweden-martina.com							Sweden&Martina	X	X	X	X	X	X			X	X	X			X	
T.F.I. System srl (IT) www.tfi-system.it							Easy Grip			X	X				X		X	X	X	X	X	-CE
Thommen Medical (CH) www.thommenmedical.com							Thommen			X	X					X	X				X	-CE
TruCrown (US) www.trucrown.com	Information not available																				X	
Vulcan Custom Dental (US) www.vulcandental.com			👍				BioHorizons	X		X	X	X	X	X	X	X	X	X	X	X	X	
Zimmer Dental (US) www.zimmerdental.com			👍					X	X					X							X	-CE
Zahn (Zirlux) (US) www.zahndental.com	Information not available																				X	

**Please note:**

- The above information is solely provided by third parties, and 3Shape assumes no responsibility for its accuracy or completeness.



**DEPARTMENT OF HEALTH & HUMAN SERVICES**

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

K151455  
3Shape A/S  
Device Trade Name: 3Shape Abutment Designer Software  
Contact Name: Hanne Nielsen

**DEFICIENCY LIST**

(b)(4)



**Page 2 - Hanne Nielsen**

(b)(4)



**Page 3 - Hanne Nielsen**

(b)(4)



**Page 4 - Hanne Nielsen**

(b)(4)



**Page 5 - Hanne Nielsen**

(b)(4)



(b)(4)



**Page 7 - Hanne Nielsen**

(b)(4)





Page 8 - Hanne Nielsen

(b)(4)





**DEPARTMENT OF HEALTH & HUMAN SERVICES**

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

K151455  
3Shape A/S  
Device Trade Name: 3Shape Abutment Designer Software  
Contact Name: Hanne Nielsen

**DEFICIENCY LIST**

(b)(4)



**Page 2 - Hanne Nielsen**

(b)(4)



**Page 3 - Hanne Nielsen**

(b)(4)



**Page 4 - Hanne Nielsen**

(b)(4)



**Page 5 - Hanne Nielsen**

(b)(4)



**Page 6 - Hanne Nielsen**

(b)(4)



**Page 7 - Hanne Nielsen**

(b)(4)





**Page 8 - Hanne Nielsen**

(b)(4)



Ms. Nielsen,

Following our March 8, 2016, phone conference to discuss the January 11, 2016, additional information you provided, the FDA has worked to develop the following interactive deficiencies which we believe as necessary to adequately describe the submission device, demonstrate the verification and validation of the software system, and ensure substantial equivalence of the submission device with respect to other previously cleared dental implant abutment and dental implant system devices. Some of the following questions are a request to provide the information discussed in our March 8 phone conference in written form for the record. The FDA respectfully requests the following information be provided by May 9, 2016.

(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)



Respectfully,

Andrew I. Steen  
Mechanical Engineer  
Senior Lead Reviewer - Endosseous Dental Implant Systems  
FDA/ODE/DAGRID/Dental Devices Branch  
10903 New Hampshire Avenue  
WO66 - G312 (NEW OFFICE)  
Silver Spring, MD 20993-0002  
phone: 301-796-6284  
fax: 301-847-8109  
[cid:image001.gif@01CD32AC.53BB7410]

---


THIS MESSAGE IS INTENDED ONLY FOR THE USE OF THE PARTY TO WHOM IT IS ADDRESSED AND MAY CONTAIN INFORMATION THAT IS PRIVILEGED, CONFIDENTIAL, AND PROTECTED FROM DISCLOSURE UNDER LAW. If you are not the addressee, or a person authorized to deliver the document to the addressee, you are hereby notified that any review disclosure, dissemination, copying, or other action based on the content of this communication is not authorized. If you have received this document in error, please immediately notify us by email or telephone. This communication is consistent with 21 CFR 10.85 (k) and constitutes an informal communication that represents my best judgment at this time but does not constitute an advisory opinion, does not necessarily represent the formal position of FDA, and does not bind or otherwise obligate or commit the agency to the views expressed.

Excellent customer service is important to us. Please take a moment to provide feedback regarding the customer service you have received:  
<https://www.research.net/s/cdrhcustomerservice?O=400&D=430&B=436&E=&S=E>



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455 – Interactive Deficiencies

(b)(4)



CONFIDENTIAL

Page 1





3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455 – Interactive Deficiencies

(b)(4)

A large, solid black rectangular redaction box covers the majority of the page's content, starting below the header and ending above the footer. The text "(b)(4)" is printed in red at the top left corner of this redacted area.

CONFIDENTIAL

Page 2



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455 – Interactive Deficiencies

(b)(4)

A large, solid black rectangular redaction box covers the majority of the page's content, starting below the header and ending above the footer.

CONFIDENTIAL

Page 3



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455 – Interactive Deficiencies

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the header and ending above the footer. The text "(b)(4)" is written in red at the top left corner of this redacted area.

CONFIDENTIAL

Page 4



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455 – Interactive Deficiencies

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the section header and ending above the footer. The text "(b)(4)" is written in red at the top left corner of this redacted area.

CONFIDENTIAL

Page 5



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455 – Interactive Deficiencies

(b)(4)

A large, solid black rectangular redaction box covers the majority of the page's content, starting below the header and ending above the footer. The text "(b)(4)" is printed in red at the top-left corner of this redacted area.


CONFIDENTIAL

Page 6



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455 - Interactive Deficiencies

(b)(4)



CONFIDENTIAL

Page 7



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455 – Interactive Deficiencies

(b)(4)

A large black rectangular redaction box covers the majority of the page content, starting below the header and ending above the footer. The text "(b)(4)" is written in red at the top left corner of this redacted area.

CONFIDENTIAL

Page 8



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455 – Interactive Deficiencies

(b)(4)

A large, solid black rectangular redaction box covers the majority of the page's content, starting below the title and ending above the footer.

CONFIDENTIAL


Page 9





3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455 - Interactive Deficiencies

(b)(4)




CONFIDENTIAL

Page 10



3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455 – Interactive Deficiencies

(b)(4)



CONFIDENTIAL

Page 11




(b)(4)

CONFIDENTIAL

Page 12

3SHAPE DENTAL SYSTEM™ SOFTWARE 510(K) K151455 – Interactive Deficiencies



(b)(4)



(b)(4)



(b)(4)



(b)(4)



(b)(4)





(b)(4)



(b)(4)

