# era HEALING AT THE SPEED OF LIGHT®

FDA/CDRH/DCC KISI816 AUG 20 2015 RECEIVED

Aug 19, 2015

Food and Drug Administration Center for Devices and Radiological Health Document Mail Center WO66 - G60910903 New Hampshire Avenue Silver Springs, MD 20993-0002

RE:

Submission Number: K151816

510k Premarket Notification **TLC 2000** 

Theralase Inc. 1945 Queen Street East Toronto, Ontario M4L 1H7 Canada

Phone: 416.699.5273 Fax: 416.699.5250

Contact: Mr. Roger Dumoulin-White

Email: rwhite@theralase.com

To Whom It May Concern:

We are responding to your eCopy Hold Letter dated 08/18/2015 regarding submission number K151816 and have enclosed a replacement eCopy flash drive.

We apologize for this inconvenience.

This eCopy is an exact duplicate of the paper copy.

Theralase Inc. is requesting market clearance for their TLC2000 Laser System for use as an adjunctive treatment for knee pain therapy in accordance with 21 CFR 807.81. Devices of this type have been classified as Class II (performance standards) medical devices according to 21 CFR 890.550. Evidence to support the substantial equivalence to a predicated device is included in this submission.

1945 Queen Street East, Toronto, Ontario M4L 1H7, Canada



The information contained herein is deemed confidential and the protection afforded to such confidentiality by application of laws and federal regulations is hereby claimed.

Your prompt response and affirmative reply to this Premarket Notification is greatly appreciated. Should you have any questions regarding this submission please do not hesitate to contact me.

Mr. Roger Dumoulin-White President and CEO Theralase Inc.

1945 Queen Street East, Toronto, Ontario M4L 1H7, Canada • Allestions?Contact 5DA/CDRH/OCE/DID at SDRH-E015TATUS.GOV@fra.Ws.efver:alases.tem Processed under FOIA Request 2016-6074. Released by



Aug 13, 2015

Food and Drug Administration Center for Devices and Radiological Health Document Mail Center WO66 – G609 10903 New Hampshire Avenue Silver Springs, MD 20993-0002

RE:

Submission Number: K151816

510k Premarket Notification TLC 2000

Theralase Inc. 1945 Queen Street East Toronto, Ontario M4L 1H7 Canada

Phone: 416.699.5273 Fax: 416.699.5250

Contact: Mr. Roger Dumoulin-White

Email: rwhite@theralase.com

To Whom It May Concern:

We are responding to your eCopy Hold Letter dated 07/06/2015 regarding submission number K151816.

We have reviewed the letter and have determined that the list of errors is due to hidden files resident on the flash drive. We apologize for this inconvenience. We have taken precautions to ensure that there are no hidden files on the included replacement flash drive.

All files are saved as pdf files which can be opened by Adobe Reader Version 9.5.

We ask that the original hard copy and eCopy submission be replaced with the hard copy and eCopy of this 510(k) submission. This eCopy is an exact duplicate of the paper copy.

# FDA CDRH DMC AUG 1 8 2015

Received

1945 Queen Street East, Toronto, Ontario M4L 1H7, Canada •416.699.5275A/CDRH/01-8666484345233ATUS.GOV@MMMstheralase.com



Theralase Inc. is requesting market clearance for their TLC2000 Laser System for use as an adjunctive treatment for knee pain therapy in accordance with 21 CFR 807.81. Devices of this type have been classified as Class II (performance standards) medical devices according to 21 CFR 890.550. Evidence to support the substantial equivalence to a predicated device is included in this submission.

The information contained herein is deemed confidential and the protection afforded to such confidentiality by application of laws and federal regulations is hereby claimed.

Your prompt response and affirmative reply to this Premarket Notification is greatly appreciated. Should you have any questions regarding this submission please do not hesitate to contact me.

Mr. Roger Dumoulin-White President and CEO Theralase Inc.

1945 Queen Street East, Toronto, Ontario M4L 1H7, Canada •416.698.5275 A/CDRH/0128662.86345233 ATUS GOV@WWW.theralase.com



June 26, 2015

Food and Drug Administration Center for Devices and Radiological Health Document Mail Center WO66 – G609 10903 New Hampshire Avenue Silver Springs, MD 20993-0002

RE: 510k Premarket Notification TLC 2000

> Theralase Inc. 1945 Queen Street East Toronto, Ontario M4L 1H7 Canada

Phone: 416.699.5273 Fax: 416.699.2250

Contact: Mr. Roger Dumoulin-White

Email: rwhite@theralase.com



JUL 0 6 2015

Received

To Whom It May Concern:

Theralase Inc. is requesting market clearance for their TLC2000 Laser System for use as an adjunctive treatment for knee pain therapy in accordance with 21 CFR 807.81. Devices of this type have been classified as Class II (performance standards) medical devices according to 21 CFR 890.550. Evidence to support the substantial equivalence to predicated devices of the same type is included in this submission.

Enclosed herein is the original hard copy and one electronic copy (eCopy) of the 510k Notification submission. This eCopy is an exact duplicate of the paper copy.

The information contained herein is deemed confidential and the protection afforded to such confidentiality by application of laws and federal regulations is hereby claimed.



A copy of the Medical Device User Fee Cover Sheet and the payment verification are enclosed for your information. These documents have been sent to the appropriate division of the Department of Health and Human Services as required.

Your prompt response and affirmative reply to this Premarket Notification is greatly appreciated. Should you have any questions regarding this submission please do not hesitate to contact me.

Mr. Roger Dumoulin-White President and CEO Theralase Inc.

1945 Queen Street East, Toronto, Ontario M4L 1H7, Canada •416.699.5273 •1.866.843.5273 •www.theralase.com Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@rda.hins.gov or 301-796-8188

# HEALING AT THE SPEED OF LIGHT®

Aug 13, 2015

Food and Drug Administration Center for Devices and Radiological Health Document Mail Center WO66 – G609 10903 New Hampshire Avenue Silver Springs, MD 20993-0002

RE:

Submission Number: K151816

510k Premarket Notification TLC 2000

Theralase Inc. 1945 Queen Street East Toronto, Ontario M4L 1H7 Canada

Phone: 416.699.5273 Fax: 416.699.5250

Contact: Mr. Roger Dumoulin-White

Email: rwhite@theralase.com

To Whom It May Concern:

We are responding to your eCopy Hold Letter dated 07/06/2015 regarding submission number K151816.

We have reviewed the letter and have determined that the list of errors is due to hidden files resident on the flash drive. We apologize for this inconvenience. We have taken precautions to ensure that there are no hidden files on the included replacement flash drive.

All files are saved as pdf files which can be opened by Adobe Reader Version 9.5.

We ask that the original hard copy and eCopy submission be replaced with the hard copy and eCopy of this 510(k) submission. This eCopy is an exact duplicate of the paper copy.

1945 Queen Street East, Toronto, Ontario M4L 1H7, Canada • 416.699.5275 • 1.866.843.5273 • WWW.theralase.com

Theralase Inc. is requesting market clearance for their TLC2000 Laser System for use as an adjunctive treatment for knee pain therapy in accordance with 21 CFR 807.81. Devices of this type have been classified as Class II (performance standards) medical devices according to 21 CFR 890.550. Evidence to support the substantial equivalence to a predicated device is included in this submission.

The information contained herein is deemed confidential and the protection afforded to such confidentiality by application of laws and federal regulations is hereby claimed.

Your prompt response and affirmative reply to this Premarket Notification is greatly appreciated. Should you have any questions regarding this submission please do not hesitate to contact me.

era

HEALING AT THE SPEED OF LIGHT®

Mr. Roger Dumoulin-White President and CEO Theralase Inc.

1945 Queen Street East, Toronto, Ontario M4L 1H7, Canada •416:699.5275<sup>DA/CDRH/OCE/DIDateDRH-E0/STATUS.GOV@Wahther301a38e.286</sup>m

#### **TABLE OF CONTENTS**


File #	Section	Description
001		Cover Letter
002		Table of Contents
003	1.1	Medical Device User Fee Cover Sheet Form FDA 3601
004	1.2	CDRH Premarket Review Submission Cover Sheet
005	1.3	Class III Summary and Certification
006	1.4	Financial Certification Statement
007	1.5	Truth and Accuracy Statement
008	1.6	510(k) Statement
009	1.7	Indications for Use Statement
010	1.8	Device Description
011	1.9	Declaration of Conformity
012	2.1	Product Description
013	2.2	Safety
014	2.3	Engineering Documentation
015	2.4	Sterilization
016	2.5	Biocompatibility and Toxicity
017	2.6	Proposed Labeling
018	2.7	Instructions for Use
019	2.8	Advertising and Promotional Materials

020	3.1	Substantial Equivalence
021	3.2	Software Verification and Validation
022	3.3	Risk Analysis
023	3.4	Predicate Device Documentation
	Appendices	
024	A	TLC 2000 User's Guide / Manual
025	В	510(k) Safety and Efficacy
026	C	MSDS
027	D	Decision Flowchart and Annex A
(b)(4)Proprie	tary Information	

#### Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017 SOPP 8407 Appendix 1 – Medical Device User Fee Cover Sheet

MEDICAL DEVICE USER FEE COVER SHEET A completed cover sheet must accompany each original application courier, please include a copy of this completed form with payment. http://www.fda.gov/oc/mdufma/coversheet.html 1. COMPANY NAME AND ADDRESS (include name, street address, city state, country, and post office code) Theralase Inc. 1945 Queen Street East Toronto, Ontario M4L 1H7, Canada 1.1 EMPLOYER IDENTIFICATION NUMBER (EIN)	2. CONTACT NAME Roger Dumoulin-White 2.1 E-MAIL ADDRESS rwhite@theralase.com 2.2 TELEPHONE NUMBER (include Area code) 416 699 5273 2.3 FACSIMILE (FAX) NUMBER (Include Area code)
address, city state, country, and post office code) Theralase Inc. 1945 Queen Street East Toronto, Ontario M4L 1H7, Canada	Roger Dumoulin-White         2.1 E-MAIL ADDRESS         rwhite@theralase.com         2.2 TELEPHONE NUMBER (include Area code)         416 699 5273         2.3 FACSIMILE (FAX) NUMBER (Include Area code)
1945 Queen Street East Toronto, Ontario M4L 1H7, Canada	<ul> <li>2.1 E-MAIL ADDRESS rwhite@theralase.com</li> <li>2.2 TELEPHONE NUMBER (include Area code) 416 699 5273</li> <li>2.3 FACSIMILE (FAX) NUMBER (Include Area code)</li> </ul>
1945 Queen Street East Toronto, Ontario M4L 1H7, Canada	rwhite@theralase.com 2.2 TELEPHONE NUMBER (include Area code) 416 699 5273 2.3 FACSIMILE (FAX) NUMBER (Include Area code)
1945 Queen Street East Toronto, Ontario M4L 1H7, Canada	<ul> <li>2.2 TELEPHONE NUMBER (include Area code)</li> <li>416 699 5273</li> <li>2.3 FACSIMILE (FAX) NUMBER (Include Area code)</li> </ul>
M4L 1H7, Canada	416 699 5273 2.3 FACSIMILE (FAX) NUMBER (Include Area code)
	2.3 FACSIMILE (FAX) NUMBER (Include Area code)
1.1 EMPLOYER IDENTIFICATION NUMBER (EIN)	
I. I EMPLOTER IDENTIFICATION NUMBER (EIN)	
	416 699 5250
descriptions at the following web site: http://www.fda.gov/oc/mdufma	ving in each column; if you are unsure, please refer to the application a 3.1 Select a center
K) Premarket notification(510(k)); except for third party	CDRH
[] 513(g) Request for Information	[] CBER
[] Biologics License Application (BLA)	3.2 Select one of the types below
[] Premarket Approval Application (DEA)	Original Application
[] Modular PMA	Supplement Types:
[] Product Development Protocol (PDP)	[] Efficacy (BLA)
[] Premarket Report (PMR)	[] Panel Track (PMA, PMR, PDP)
[] Annual Fee for Periodic Reporting (APR)	[] Real-Time (PMA, PMR, PDP)
[] 30-Day Notice	[] 180-day (PMA, PMR, PDP)
<ul> <li>4. ARE YOU A SMALL BUSINESS? (See the instructions for more [] YES, I meet the small business criteria and have submitted the r qualifying documents to FDA</li> <li>4.1 If Yes, please enter your Small Business Decision Number:</li> </ul>	
5. FDA WILL NOT ACCEPT YOUR SUBMISSION IF YOUR COMP THAT IS DUE TO FDA. HAS YOUR COMPANY PAID ALL ESTABL [4] YES (All of our establishments have registered and paid the fee, 30 days of FDA's approval/clearance of this device.)	ANY HAS NOT PAID AN ESTABLISHMENT REGISTRATION FEE LISHMENT REGISTRATION FEES THAT ARE DUE TO FDA? , or this is our first device, and we will register and pay the fee within paid all fees due to FDA. This submission will not be processed; see
http://www.fda.gov/cdrh/mdufma for additional information)	
APPLICABLE EXCEPTION. No	THE FOLLOWING USER FEE EXCEPTIONS? IF SO, CHECK THE
[] This application is the first PMA submitted by a qualified small business, including any affiliates	[] The sole purpose of the application is to support conditions of use for a pediatric population
[] This biologics application is submitted under section 351 of the F Health Service Act for a product licensed for further manufacturing u only	
7. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION F PEDIATRIC POPULATION THAT NOW PROPOSES CONDITION subject to the fee that applies for an original premarket approval app	OF USE FOR ANY ADULT POPULATION? (If so, the application is

Payment was not made using the Medical Device User Fee Cover Sheet Form 3601. Consequently, there is no MD number on the Medical Device User Fee Cover Sheet. However, the proof of payment information auto-reply is provided below.



# Please wait...

If this message is not eventually replaced by the proper contents of the document, your PDF viewer may not be able to display this type of document.

You can upgrade to the latest version of Adobe Reader for Windows®, Mac, or Linux® by visiting http://www.adobe.com/products/acrobat/readstep2.html.

For more assistance with Adobe Reader visit http://www.adobe.com/support/products/ acrreader.html.


Windows is either a registered trademark or a trademark of Microsoft Corporation in the United States and/or other countries. Mac is a trademark of Apple Inc., registered in the United States and other countries. Linux is the registered trademark of Linus Torvalds in the U.S. and other countries.

#### 1.3 Class III Summary and Certification

The TLC-2000 is Class II, Performance Standards, 21CFR Part 890.5500 – Infrared Lamp, Non-Heating. Consequently this section does not apply to the TLC-2000 product.

#### 1.4 Financial Certification Statement

This section is not applicable to this submission as there were no clinical studies commissioned.



## Premarket Notification Truthful And Accurate Statement [As Required by 21 CFR 807.87(k)]

I certify that, in my capacity as President and CEO of

Theralase Inc., I believe to the best of my knowledge, that all data

and information submitted in the premarket notification are truthful and

accurate and that no material fact has been omitted.

(Signature)
Roger Dumoulin-White
(Typed Name) August 15, 2015
(Date)
K151816

ieral

HEALING AT THE SPEED OF LIGHT®

\*(Premarket Notification [510(k)] Number)

\*For a new submission, leave the 510(k) number blank.

Must be signed by a responsible person of the firm required to submit the premarket notification [e.g., not a consultant for the 510(k) submitter].

1945 Queen Street East, Toronto, Ontario M4L 1H7, Canada • 49.69.699.992975 DA/CDRH/OCE/018-85478-521954TUS.GOV@ff4-W9:8916F314388.268m

# Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017 HEALING AT THE SPEED OF LIGHT®

### Premarket Notification 510(k) Statement

(As Required By 21 CFR 807.93)

I certify that, in my capacity as President and CEO of Theralase Inc., I will make available all information included in this premarket notification on safety and effectiveness within 30 days of request by any person if the device described in the premarket notification submission is determined to be substantially equivalent. The information I agree to make available will be a duplicate of the premarket notification submission, including any adverse safety and effectiveness information, but excluding all patient identifiers, and trade secret and confidential commercial information, as defined in 21 CFR 20.61.

(Signature of Certifier)

Roger Dumoulin-White (Typed Name)

Juguet 13, 2015 (Date)

K151816 \*(Premarket Notification [510(k)] Number)

\*For a new submission, leave the 510(k) number blank.

1945 Queen Street East, Toronto, Ontario M4L 1H7, Canada • 4Questions 9 Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov.or 301-796-8188 DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

#### Indications for Use

510(k) Number *(if known)* Pending

Device Name TLC-2000 Therapeutic Medical Laser System

Indications for Use (Describe)

The Theralase TLC-2000 Therapeutic Medical Laser System is indicated for adjunctive use in the temporary relief of pain associated with knee disorders. An average power setting of 40 mW is used with a treatment time of 1 minute per location. There are seven locations around the knee, three each lateral and medial. The seventh location is on the posterior of the knee at midline of the popliteal fossa.

Type of Use (Select one or both, as applicable)	
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)

#### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov* 

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Form Approved: OMB No. 0910-0120 Expiration Date: January 31, 2017 See PRA Statement below.

EF

#### **1.8 Product Description**

The Theralase<sup>TM</sup> TLC-2000 Therapeutic Medical Laser System is a Class 3B medical laser system and is indicated for use for adjunctive use in the temporary relief of pain associated with knee disorders. The Theralase<sup>TM</sup> TLC-2000 Therapeutic Medical Laser System is comprised of the following components:

b)(4)Proprietary Information

The power pack only supplies power to the probe(s) and incorporates the necessary safety measures for laser equipment the emergency stop, key switch and remote shunt. The probes are controlled by Theralase software resident on the tablet computer. The practitioner selects the knee pain diagnosis and the Theralase software transmits the correct power level and time settings for the knee pain protocol to the probe(s). Communications between the tablet and the probe(s) is via a wireless Blue Tooth link. (b)(4)Proprietary Information

The TLC-2000 laser probes are used in direct contact with tissue in order to inject photons of light (b)(4)Proprietary Information non-invasively and drug-free into tissue.

The external power supply is a medically approved power supply. The input voltage range is 100 240 VAC with a frequency range of 47 73 Hz. The power supply connects to the AC source with a medically approved ac cord.

The system is supplies with the Operations Manual, Pain Management Manual, safety glasses and a carrying case.

The Operations Manual is attached as Appendix A to this submission.

#### Section 1.9 Declaration of Conformity and Summary Reports (FORM FDA 3654)

The TLC2000 Laser System is designed to comply with the Federal Food and Drug Administration's Code of Federal Regulations 21 CFR 1010 and 21 CFR 1040 and other applicable federal and state performance standards and regulations. The Company declares that the technical information provided within this submission is in conformance with the following consensus standards:

(b)(4)Proprietary Information

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017 Form Approved: OMB No. 0910-0120; Expiration Date: 1/31/2017

Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)				
This report and the Summary Report Table are to be comp ences a national or international standard. A separate report				
TYPE OF 510(K) SUBMISSION	Abbreviated			
STANDARD TITLE ' IEC 60825-1: Safety of Laser Products, Part 1: Equipment Class	sification and Requirements			
Please answer the following questions		Yes	No	
Is this standard recognized by FDA <sup>2</sup> ?		$\boxtimes$		
FDA Recognition number <sup>3</sup>		# <u>12-273</u>		
Was a third party laboratory responsible for testing conform in the 510(k)?		$\boxtimes$		
Is a summary report <sup>4</sup> describing the extent of conformance 510(k)? If no, complete a summary report table.		$\boxtimes$		
Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device?				
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).				
Does this standard include more than one option or selectic If yes, report options selected in the summary report table.	on of tests?		$\boxtimes$	
Were there any deviations or adaptations made in the use of If yes, were deviations in accordance with the FDA supplem			$\square$	
Were deviations or adaptations made beyond what is speci If yes, report these deviations or adaptations in the summar			$\boxtimes$	
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			$\boxtimes$	
Is there an FDA guidance <sup>6</sup> that is associated with this stand If yes, was the guidance document followed in preparation of Title of guidance:	of this 510k?			
<ol> <li><sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]</li> <li><sup>2</sup> Authority [21 U.S.C. 360d], http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm</li> <li><sup>3</sup> http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm</li> <li><sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the</li> </ol>	address of the test laboratory or certification body inv assessment to this standard. The summary report inc all standards utilized during the development of the d <sup>5</sup> The supplemental information sheet (SIS) is addition. is necessary before FDA recognizes the standard. Fo www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStand <sup>6</sup> The online search for CDRH Guidance Documents c http://www.fda.gov/MedicalDevices/DeviceRegulation GuidanceDocuments/default.htm	cludes inforr evice. al informatic bund at http ards/search an be found	mation on on which :// n.cfm I at	

EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE					
STANDARD TITLE IEC 60825-1: Safety	of Laser Products, Part 1: Equipment C	lassification and Requirements			
	CONFORMANCE WIT	TH STANDARD SECTIONS*			
SECTION NUMBER	SECTION TITLE		CONFORMANCE?		
			Yes No N/A		
TYPE OF DEVIATION OF	OPTION SELECTED *				
DESCRIPTION					
DESCRIPTION					
JUSTIFICATION					
SECTION NUMBER	SECTION TITLE		CONFORMANCE?		
			Yes No N/A		
TYPE OF DEVIATION OF	COPTION SELECTED				
DESCRIPTION					
JUSTIFICATION					
SECTION NUMBER	SECTION TITLE		CONFORMANCE?		
TYPE OF DEVIATION OF	OPTION SELECTED *		Yes No N/A		
DESCRIPTION					
JUSTIFICATION					
explanation is neede described and adequ selected when follow	d under "justification." Some standards in ately justified as appropriate for the subje	whether conformance is met. If a section include options, so similar to deviations, th act device. Explanation of all deviations of deviation or option selected," "description	e option chosen needs to be r description of options		
		ne standard, a deviation brought out by th he device, or any adaptation of a section.			
		ents of the Paperwork Reduction Act of 1995			
		I TO THE PRA STAFF EMAIL ADDRES			
instructions, search information. Send	existing data sources, gather and main	ated to average 1 hour per response, in ntain the data needed and complete and ate or any other aspect of this informa-	d review the collection of		
Departm Food an Office o	ent of Health and Human Services d Drug Administration f Chief Information Officer	"An agency may not con a person is not require collection of information	ed to respond to, a 1 unless it displays a		
	ork Reduction Act (PRA) Staff ff@fda.hhs.gov	currently valid OMB			
FORM FDA 3654 (4/14	)	Page 2 of 2	Add page		

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017 Form Approved: OMB No. 0910-0120; Expiration Date: 1/31/2017

Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)					
This report and the Summary Report Table are to be comp ences a national or international standard. A separate repor					
TYPE OF 510(K) SUBMISSION	Abbreviated				
STANDARD TITLE <sup>1</sup> IEC 60601-1:2012 Medical Electrical Equipment, Part 1: Genera	al Requirements for Basic Safety and Essential	Perform	ance		
Please answer the following questions		Yes	No		
Is this standard recognized by FDA <sup>2</sup> ?		$\boxtimes$			
FDA Recognition number <sup>3</sup>		# <u></u> 5-76			
Was a third party laboratory responsible for testing conform in the 510(k)?		$\boxtimes$			
Is a summary report <sup>4</sup> describing the extent of conformance 510(k)? If no, complete a summary report table.					
Does the test data for this device demonstrate conformity to pertains to this device?		$\boxtimes$			
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).		$\boxtimes$			
Does this standard include more than one option or selection If yes, report options selected in the summary report table.	n of tests?		$\boxtimes$		
Were there any deviations or adaptations made in the use of If yes, were deviations in accordance with the FDA supplem					
Were deviations or adaptations made beyond what is speci If yes, report these deviations or adaptations in the summar			$\boxtimes$		
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			$\boxtimes$		
Is there an FDA guidance <sup>6</sup> that is associated with this stand If yes, was the guidance document followed in preparation of Title of guidance:					
<ul> <li><sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]</li> <li><sup>2</sup> Authority [21 U.S.C. 360d], http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm</li> <li><sup>3</sup> http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm</li> <li><sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and</li> </ul>	address of the test laboratory or certification body inv assessment to this standard. The summary report inc all standards utilized during the development of the d <sup>5</sup> The supplemental information sheet (SIS) is addition is necessary before FDA recognizes the standard. Fo www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStand <sup>6</sup> The online search for CDRH Guidance Documents c http://www.fda.gov/MedicalDevices/DeviceRegulation GuidanceDocuments/default.htm	cludes info evice. al informat bund at http ards/searc an be foun	rmation on ion which p:// h.cfm d at		

		ANDARD CONFORMANCE RY REPORT TABLE	
STANDARD TITLE IEC 60601-1:2012	Medical Electrical Equipment, Part 1:	General Requirements for Basic Safety and	Essential Performance
	CONFORMANCE	WITH STANDARD SECTIONS*	
SECTION NUMBER	SECTION TITLE		CONFORMANCE?
			Yes No N/A
TYPE OF DEVIATION (	DR OPTION SELECTED *		1
DESCRIPTION			
JUSTIFICATION			
SECTION NUMBER	SECTION TITLE		CONFORMANCE?
	ž		
TYPE OF DEVIATION (	DR OPTION SELECTED *		
DESCRIPTION			
JUSTIFICATION			
SECTION NUMBER	SECTION TITLE		CONFORMANCE?
			Yes No N/A
TYPE OF DEVIATION C	R OPTION SELECTED *		
DESCRIPTION			
JUSTIFICATION			
explanation is need described and adeq selected when follow	ed under "justification." Some standards uately justified as appropriate for the su	ate whether conformance is met. If a section s include options, so similar to deviations, th bject device. Explanation of all deviations o of deviation or option selected," "description	e option chosen needs to be r description of options
* Types of deviations information sheet (S	can include an exclusion of a section in SIS), a deviation to adapt the standard to	the standard, a deviation brought out by th the device, or any adaptation of a section.	e FDA supplemental
	This section applies only to require	ements of the Paperwork Reduction Act of 1995.	
*D(	) NOT SEND YOUR COMPLETED FOI	RM TO THE PRA STAFF EMAIL ADDRES	S BELOW.*
instructions, searce information. Send suggestions for rea	h existing data sources, gather and m l comments regarding this burden esti ducing this burden, to:	imated to average 1 hour per response, inc aintain the data needed and complete and imate or any other aspect of this informa	l review the collection of
Food an Office o Paperw	nent of Health and Human Services ad Drug Administration of Chief Information Officer ork Reduction Act (PRA) Staff aff@fda.hhs.gov	"An agency may not cond a person is not require collection of information currently valid OMB o	d to respond to, a unless it displays a

Carps Approved	OMD No.	0010 0100.	Eunisotion	Data	1/01/0017
Form Approved	. UND NO.	. 0910-0120,	Expiration	Dale.	1/31/2017

Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)				
This report and the Summary Report Table are to be completed by the applicant when submitting a ences a national or international standard. A separate report is required for each standard referenced i				
TYPE OF 510(K) SUBMISSION				
STANDARD TITLE <sup>1</sup> IEC 60601-1-2:2007 Medical Electrical Equipment, General Requirement for Safety. Electromagnetic Compati	bility			
Please answer the following questions	Yes	No		
Is this standard recognized by FDA <sup>2</sup> ?	$\boxtimes$			
FDA Recognition number <sup>3</sup>	<u>4</u> 19-1	0		
Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)?	$\boxtimes$			
Is a summary report <sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? If no, complete a summary report table.	$\boxtimes$			
Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device?	$\boxtimes$			
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).	$\boxtimes$			
Does this standard include more than one option or selection of tests? If yes, report options selected in the summary report table.		$\boxtimes$		
Were there any deviations or adaptations made in the use of the standard? If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) <sup>5</sup> ?		$\square$		
Were deviations or adaptations made beyond what is specified in the FDA SIS? If yes, report these deviations or adaptations in the summary report table.		$\boxtimes$		
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.		$\boxtimes$		
Is there an FDA guidance <sup>6</sup> that is associated with this standard? If yes, was the guidance document followed in preparation of this 510k? Title of guidance:				
<ul> <li><sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]</li> <li><sup>2</sup> Authority [21 U.S.C. 360d], http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm</li> <li><sup>3</sup> http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm</li> <li><sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods): choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and</li> <li><sup>a</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods): choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and</li> </ul>	ludes infor evice. al informati ound at http ards/searc an be foun	mation on on which b:// h.cfm d at		

PSC Publishing Services (301) 443-6740 EF

		NDARD CONFORMANCE Y REPORT TABLE			
STANDARD TITLE IEC 60601-1-2:2007	Medical Electrical Equipment, Genera	l Requirement for Safety. Electromagneti	c Compatib	ility	
	CONFORMANCE W	TH STANDARD SECTIONS*			
SECTION NUMBER	SECTION TITLE		CONFORM	ANCE?	
			Yes	No	□ N/A
TYPE OF DEVIATION OR	OPTION SELECTED *		]		
DESCRIPTION					
JUSTIFICATION					
SECTION NUMBER	SECTION TITLE		CONFORM	ANCE2	
SECTION NOMBER			☐ Yes		□ N/A
TYPE OF DEVIATION OR	OPTION SELECTED *				
DESCRIPTION					
JUSTIFICATION					
SECTION NUMBER	SECTION TITLE			IANCE?	∏ N/A
TYPE OF DEVIATION OR	OPTION SELECTED *				
DESCRIPTION					
JUSTIFICATION					
explanation is needed described and adequ selected when followi	d under "justification." Some standards i ately justified as appropriate for the sub	e whether conformance is met. If a sectior nclude options, so similar to deviations, th ject device. Explanation of all deviations of f deviation or option selected," "description	ne option cho or descriptio	osen nee n of optic	eds to be ons
		he standard, a deviation brought out by the device, or any adaptation of a section.		olementa	41
		nents of the Paperwork Reduction Act of 1995			
		M TO THE PRA STAFF EMAIL ADDRES			
instructions, search information. Send	existing data sources, gather and ma	nated to average 1 hour per response, in intain the data needed and complete an nate or any other aspect of this inform	d review th	ne collec	tion of
Food and Office of Paperwo	ent of Health and Human Services d Drug Administration f Chief Information Officer rk Reduction Act (PRA) Staff f@fda.hhs.gov	"An agency may not con a person is not require collection of information currently valid OMB	ed to respond n unless it dis	d to, a splays a	

	Form Approved: OMB No. 0910-0120; Exp	iration Da	te: 1/31/201		
Department of Health and Human Services Food and Drug Administration STANDARDS DATA REPORT FOR 510(k)s (To be filled in by applicant)					
This report and the Summary Report Table are to be comp ences a national or international standard. A separate repor	pleted by the applicant when submitting a rt is required for each standard referenced	510(k) t in the 51	hat refer- 10(k).		
TYPE OF 510(K) SUBMISSION					
Traditional Special	Abbreviated				
STANDARD TITLE ' IEC 60601-2-22:2007 Medical Electrical Equipment, Part 2, Par Laser Equipment.	ticular Requirements for the Safety of diagnost	ic and T	herapeutic		
Please answer the following questions		Yes	No		
Is this standard recognized by FDA <sup>2</sup> ?		$\boxtimes$			
FDA Recognition number <sup>3</sup>		# <u>12-208</u>			
Was a third party laboratory responsible for testing conform in the 510(k)?	ity of the device to this standard identified	$\boxtimes$			
Is a summary report <sup>4</sup> describing the extent of conformance 510(k)? If no, complete a summary report table.					
Does the test data for this device demonstrate conformity to pertains to this device?	the requirements of this standard as it	$\boxtimes$			
Does this standard include acceptance criteria? If no, include the results of testing in the 510(k).		$\boxtimes$			
Does this standard include more than one option or selectio If yes, report options selected in the summary report table.	n of tests?		$\boxtimes$		
Were there any deviations or adaptations made in the use of If yes, were deviations in accordance with the FDA supplem			$\square$		
Were deviations or adaptations made beyond what is specif If yes, report these deviations or adaptations in the summar			$\boxtimes$		
Were there any exclusions from the standard? If yes, report these exclusions in the summary report table.			$\boxtimes$		
Is there an FDA guidance <sup>6</sup> that is associated with this stand If yes, was the guidance document followed in preparation of Title of guidance:					
Title of guidance:					
<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]	address of the test laboratory or certification body inv assessment to this standard. The summary report inc	ludes infor			
<sup>2</sup> Authority [21 U.S.C. 360d], http://www.fda.gov/MedicalDevices/ DeviceRegulationandGuidance/Standards/default.htm	all standards utilized during the development of the de <sup>5</sup> The supplemental information sheet (SIS) is additional is necessary before FDA recognizes the standard. For www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStanda				
<sup>3</sup> http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm					
<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and	Image review (for example, alternative test methods): choices made of the online search for CDRH Guidance Documents can be for http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidanceDocuments can be for http://www.fda.gov/MedicalDevi				

4/14) Page 1 of 2 Psc P Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188

PSC Publishing Services (301) 443-6740 EF

EXTENT OF STANDARD CONFORMANCE SUMMARY REPORT TABLE						
STANDARD TITLE         IEC 60601-2-22:2007       Medical Electrical Equipment, Part 2, Particular Requirements for the Safety of diagnostic and Therapeutic         Laser Equipment.						
	CONFORMANCE WITH STANDARD SECTIONS*					
SECTION NUMBER	SECTION TITLE	CONFORM	ANCE?			
		Yes	No	N/A		
TYPE OF DEVIATION OF	OPTION SELECTED *					
DESCRIPTION						
JUSTIFICATION						
SECTION NUMBER	SECTION TITLE		IANCE? □No			
TYPE OF DEVIATION OF	COPTION SELECTED *					
DESCRIPTION						
JUSTIFICATION						
SECTION NUMBER	SECTION TITLE	CONFORM	ANCE?			
TYPE OF DEVIATION OR	OPTION SELECTED *	Yes	No	N/A		
DESCRIPTION						
JUSTIFICATION						
* For completeness list all sections of the standard and indicate whether conformance is met. If a section is not applicable (N/A) an explanation is needed under "justification." Some standards include options, so similar to deviations, the option chosen needs to be described and adequately justified as appropriate for the subject device. Explanation of all deviations or description of options selected when following a standard is required under "type of deviation or option selected," "description" and "justification" on the report. More than one page may be necessary.						
* Types of deviations can include an exclusion of a section in the standard, a deviation brought out by the FDA supplemental information sheet (SIS), a deviation to adapt the standard to the device, or any adaptation of a section.						
	This section applies only to requirements of the Paperwork Reduction Act of 199	95.				
*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.*						
The burden time for this collection of information is estimated to average 1 hour per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to: Department of Health and Human Services						
Department of Health and Human Services"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number."Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number."						

#### 2.1 Device Description

#### a) THE TLC-2000 SYSTEM

a) THE ILC-2000 STSTEM	
(b)(4)Proprietary Information	

### System Component Description:

The Probe:

b)(4)Proprietary Information



The Tablet:




- 1 Power On
- 2 Charging Port
- 3 USB Port
- 4 Release Button
- 5 Touch Pad

The tablet is a windows based computer running Windows 8.1 operating system. The tablet has a Blue Tooth 4.0 radio and WiFi. The tablet has a 10.1" touch screen and a detachable keyboard. The Theralase application software is loaded on the tablet when supplied with the system.



#### **The Power Pack:**



- 1 Remote Interlock
- 2 Emergency Stop
- 3 Key Switch
- 4 Power Input
- 5 Probe Power Ports
- 6 Power In / Power Out LED indicators

The power pack has the electronics which provide the operating voltages for the probe. The power pack contains the necessary safeties in the form of an emergency stop switch, a remote interlock and a key operated on/off switch. There are four ports available to connect up to four probes.

(b)(4)Proprietary Information

#### The Power Supply:

The power supply is a medically approved universal, 100 240 VAC input voltage rating at 47/63 Hz, power supply. The output is 24 volts DC for the power pack. A medically approved ac power cord is included.



#### b) ACCESSORY COMPONENT DESCRIPTIONS

#### Safety eyewear:

Supplied with the system are two sets of standard, TLC-901 laser safety glasses one pair for the clinician and one pair for the patient. These laser safety glasses are provided for viewing diffused laser light only and are not for direct viewing of the laser output.





### c) Specifications

(b)(4)Proprietary Information	

#### d) Intended Use of the TLC-2000

The TLC-2000 therapeutic medical laser probes are used for adjunctive use in the temporary relief of pain associated with knee disorders. The technology of Low Level Laser Therapy (LLLT) is designed as a non-invasive and non-toxic therapeutic treatment. The TLC-2000 laser probes are used in direct contact with tissue in order to inject photons of light (b)(4)Proprietary Information non-invasively and drug-free into tissue.

#### e) Principles of operation

(b)(4)Proprietary Information

b)(4)Proprietary Information

(b)(4)Proprietary Information

(b)/A) Dreprietory ( Information	
(b)(4)Proprietary Information	

#### 2.2 Safety

There are no known significant safety or adverse effects reported in the literature from the use of Low Level Laser Therapy (LLLT). During the past ten years the predicate device has been on the market, there hav been no significant adverse effects reported from the use of the predicate device, the Theralase TLC-1000 Therapeutic Medical Laser System. As the new device, the Theralase TLC-2000, is comprised of the same materials, the same wavelengths and same physical design and number of laser diodes, no significant adverse effects are expected from the use of the new device.

However, following is a list of potential incidents that may occur during laser treatment:

- 1. Eye pain, retina or other ophthalmic damage may be noted if the laser beam is pointed directly at or into the eye for prolonged periods of time. Laser safety glasses, absorbent for the 905 nm wavelength, must be worn during all treatment sessions by both the subject and the health care professional(s).
- 2. Transient redness at the treatment site
- 3. Swelling at or near the area of laser irradiation
- 4. Tingling or itching
- 5. Temporary paresthesia

The FDA has classified bio-stimulation laser for pain attenuation as a non-significant risk device (HSS Publican 92-4159 IDE Manual, page 34). There are no greater risks associated with the use of LLLT than with any other treatment modality generally accepted for the treatment of knee disorders. The potential benefit of pain relief, decreased inflammation, increased function and improved quality of life may out-weigh any temporary adverse effect that may be observed from the use of the non-invasive Low Level Laser Therapy.

Over the past decades, an abundance of research studies, both *in vitro* and *in vivo*, have been conducted in laboratory and clinical environments. The vast majority of these studies report clear evidence that low level laser irradiation alters the cellular processes in a non-thermal and wavelength specific manner. Review of this literature does not reveal significant risks to study subjects and no unusual or adverse effects were reported during the clinical research conducted in which the predicate device was used.

The TLC-2000 is manufactured in an ISO 9001:20085 and ISO 13485:2003 compliant manufacturing facility in compliance with International ISO 9001 quality standard programs and the FDA Quality Systems Regulations (cGMP) 21 CFR Part 820.

The TLC-2000 sysem is tested and verified to meet all output technical specifications via 100% final quality control procedures.

The TLC-2000 System has been certified to the following safety standards:

# 2.3 Engineering Documentation

The Theralase TLC-2000 Therapeutic Medical Laser System is comprised of two parts that are manufactured by Theralase the laser probe and the power pack. (b)

(4)Proprietar y Information

10/41/Proprietary Information - Design         A         A         B <t< th=""><th>1</th><th>2</th><th>3</th><th>4</th><th>5</th><th>6</th></t<>	1	2	3	4	5	6
$ \begin{array}{ c c c c } & & & & & & & & & & & & & & & & & & &$	(b)(4)Proprietary Information	on - Design				
$ \begin{array}{ c c c } & & & & & & & & & & & & & & & & & & &$						
$ \begin{array}{ c c c c } & & & & & & & & & & & & & & & & & & &$						
$ \begin{array}{ c c c } & & & & & & & & & & & & & & & & & & &$						
$ \begin{array}{ c c c c } & & & & & & & & & & & & & & & & & & &$						
$ \begin{array}{ c c c c } & & & & & & & & & & & & & & & & & & &$						
$ \begin{array}{ c c c c } & & & & & & & & & & & & & & & & & & &$						
$ \begin{array}{ c c c c } \hline \begin{tabular}{ c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$						
$ \begin{array}{ c c c c } \hline \begin{tabular}{ c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$						
$ \begin{array}{ c c c c } \hline \begin{tabular}{ c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$						
$ \begin{array}{ c c c c } \hline \begin{tabular}{ c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$						
$ \begin{array}{ c c c } \hline \begin{tabular}{ c c } \hline \begi$						
$ \begin{array}{ c c c } \hline \begin{tabular}{ c c } \hline \begi$						
$ \begin{array}{ c c c } \hline \begin{tabular}{ c c } \hline \begi$						
$ \begin{array}{ c c c c } \hline \begin{tabular}{ c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c } \hline \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$						
$ \begin{array}{ c c c } \hline \begin{tabular}{ c c } \hline \begi$	в					
$ \begin{array}{ c c c } \hline \begin{tabular}{ c c } \hline \begi$						
$ \begin{array}{ c c c } \hline \begin{tabular}{ c c } \hline \begi$						
$ \begin{tabular}{ c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$						
$ \begin{tabular}{ c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$						
$ \begin{tabular}{ c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$						
$ \begin{tabular}{ c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$						
$ \begin{tabular}{ c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$						
$ \begin{tabular}{ c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$						
$ \begin{tabular}{ c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$						
$ \begin{array}{ c c c } \hline \begin{tabular}{ c c } \hline \begi$						
$ \begin{array}{ c c c } \hline \begin{tabular}{ c c } \hline \begi$	c					
$ \begin{array}{ c c c } \hline \begin{tabular}{ c c } \hline \begi$						
$ \begin{tabular}{ c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$						
$ \begin{array}{ c c c } \hline \begin{tabular}{ c c } \hline \begi$						
$ \begin{array}{ c c c } \hline \begin{tabular}{ c c } \hline \begi$						
$ \begin{tabular}{ c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$						
$ \begin{tabular}{ c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$						
$ \begin{array}{ c c c } \hline \begin{tabular}{ c c } \hline \begi$						
$ \begin{tabular}{ c c c c c } \hline \begin{tabular}{ c c c c c c c } \hline \begin{tabular}{ c c c c c c c c c c c c c c c c c c c$						
D Update						
Update:     Update:     CHECKED     Image:						9/2014 Schematic, TLC2000
Image:	D			Update:		
SIZE     NUMBER     REV       B     500-0002.     R4						Controller
SHEET 1 OF 10					Q.A	
1         2         3         4         5         6						B 500-0002. R4
1 2 3 4 5 6						SHEET 1 OF 10
	1	2	3	4	5	

#### 2.3.1\_Controller Board Schematic

1	2	3	4	5	6
(b)(4)Proprietary Information	- Design				
			Update: 20/02/2015 Update:	DRAWN JACK LIN 12/09/20	
			Update:	CHECKED	Schematic, TLC2000
D				ENG APPR.	Controller <sup>D</sup>
				Q.A	SIZE NUMBER REV
					B 500-0002. R4
1	2	3	4	5	SHEET 2 OF 10 6
	-				

	1	2	3	4	5	6
	b)(4)Proprietary Information	n - Design	• •			·
Α						
в						
C						
				Update: 20/02/2015	NAME D	A 117
				Update: 20/02/2015 Update:		
				Update:	CHECKED 120	Schematic, TLC2000
D				Gpanne.	 ENG APPR.	Controller D
					QA	
						SIZE NUMBER REV
						B 500-0002. R4
		r				SHEET 3 OF 10
	1	2	3	4	5	6

Γ		ب م ا				D	_				]
		Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017	•				1000	R4 KEV	9 80	5	
				0000	0007				SHEET 4	1	
	9			L L L	лс, 11.C	ler		500-0002.		°	
					Schemalic, 1LC2000	Controller	41000				
			DATE								
			NAME	IACK LIN							
	s			DRAWN	CHECKED	ENG APPR.	δA			5	
						Τ				1	
-	_										
			20/02/2015						1		
	4		Update: 20		Update:					4	
									1		
		nally									
-	_	Left blank intentionally								-	
		k int									
		ft blar									
	3	T								3	
	7									2	
	_										
	-									-	
		Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-79	96-	818	38						
		ب کې				D					

	1	2	3	4		5	6
	(b)(4)Proprietary Informatio	n - Dəsign					
Λ							
в							
	-						
с							
				Update: 20/02/201	5		DATE
				Update:	-	I I I	Schematic, TLC2000
D				Update:		CHECKED ENG APPR.	Controller D
						Q.A	
							SIZE         NUMBER         REV           B         500-0002.         R4
	1	2	3	4		5	SHEET 5 OF 10 6
	-	-					

1	2	3	4	5	6
(b)(4)Proprietary Information	on - Design				
A					
в					
c					
			Update: 20/02/2015	NAME	DATE
			Update: 20/02/2015		
			Update:	CHECKED	Schematic, The2000
D				ENG APPR.	Controller D
				Q.A	SIZE NUMBER REV
					B 500-0002. R4
					SHEET 6 OF 10
1	2	3	4	5	6

	1	2	3	4	5	6
~	(b)(4)Proprietary Informatio	n - Design				A
в						в
c						c
D				Update: 20/02/2015 Update: Update:		DATE 19/2014 Schematic, TLC2000 Controller D SIZE NUMBER REV
	1	2	3	4	5	B 500-0002. R4 <u>SHEET 7 OF 10</u> 6

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188

1	2	3	4		5	6
(b)(4)Proprietary Information	n - Design					
			Update:	I	DRAWN JACK LIN 13/11/2014	Schematic, TLC2000
D			Update:		CHECKED	
					ENG APPR.	Controller
					Q.A	
						SIZE NUMBER REV
						B 500-0002. R4
						SHEET 8 OF 10
1	2	3	4	· · ·	5	6
*	2		7		-	÷


1 2	3	4	5	6
(b)(4)Proprietary Information - Design	-	-	-	-
(b)(1) repristally merinament beergin				
Α				
В				
с				
D				

1	2	3	4	5	6
(b)(4)Proprietary Information	ı - Design				
-					
			Update: 20/02/2015	NAME D	ATE
			Update:		Schematic, TLC2000
			Update:	CHECKED	Schematic, TLC2000
D				ENG APPR.	Controller D
				QA	SIZE NUMBER REV
					B 500-0002. R4
					SHEET 10 OF 10
1	2	3	4	5	6

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188

2.3.2\_Controller Board PCB


Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017



CAM350 V 10.1.0 : Thu Feb 19 14:08:52 2015 - (Untitled) : ASP-600-0002-R4.pho

2.3.2\_Controller Board PCB

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017



CAM350 V 10.1.0 : Thu Feb 19 14:08:52 2015 - (Untitled) : ASS-600-0002-R4.pho

(b)(4)Proprietary Information - Design				
		TY SYN	-	Γ
		434 + 6 ×	YES	╀
		5 🖸	YES	+
	18	 29 🗇	YES	T
	F	 6 🖂 42 + <sup>4</sup>	-	╀
THERALASE TECHNOLOGIES INC	23	 42 + <sup>4</sup> 14 + <sup>8</sup>	YES	+

THERALASE TECHNOLOGIES INC				
DRILL GRAPHICS				
PROD: TLC2000 CONTROLLER REV4 DATE: Feb 19 2015 No. of PCB LAYERS: 6				
NO. OF FUE LATENS. 0				

	434	+	YES	+/003
	6	×	YES	+/003
	5	Ð	YES	+/003 TO BE FILLED
18	29	¢	YES	+/003
27.56	6	Χ	YES	+/003
28	42	+	YES	+/003
31.5	14	۳	YES	+/003
35	1	¢	YES	+/003
40	2	₽	YES	+/003
43	5	ŧ	YES	+/003
47.24	2	+	YES	+/003
125.98	4	+e	YES	+/003
35.43	2	4	NO	+/003

TOL

Refer to : PCB-600-0002-R4.DOC for Fabrication Notes

b)(4)Proprietary Information - Design	

THERALASE TECHNOLOGIES INC
LAYER 1 – PRIMARY SIDE
PROD: TLC2000 CONTROLLER REV4 DATE: Feb 19 2015 No. of PCB LAYERS: 6

CAM350 V 10.1.0 : Thu Balestibes?120nfact:FD2/CDRHJOEE/DIDat CDRHEFOISTATUS-COV@fda.6hs-gov/or/300-796-8188

(b)(4)Proprietary Information - Design	

THERALASE TECHNOLOGIES INC	
LAYER 2 - GND	
PROD: TLC2000 CONTROLLER REV4 DATE: Feb 19 2015 No. of PCB LAYERS: 6	



CAM350 V 10.1.0 : Questions? Contact #DA/CDRD/QCE/DID-at CDRH-E01STATUS.GDV@fdacture.gov.cor 3R147.96+8188



THERALASE TECHNOLOGIES INC
LAYER 3 - SIG1
PROD: TLC2000 CONTROLLER REV4 DATE: Feb 19 2015 No. of PCB LAYERS: 6



CAM350 V 10.1.0 : Questions Contact FDA/GDBH/OCEADD at CORN-FOISTATUS GDY @ Ida hts gov of 301479678188



THERALASE TECHNOLOGIES INC
LAYER 4 – SIG2
PROD: TLC2000 CONTROLLER REV4 DATE: Feb 19 2015 No. of PCB LAYERS: 6

CAM350 V 10.1.0 : PHESTIPPES CPOTAGE FDAGEDEN/OCE/BID\_at CORH-FOISTATUS. GOV @ 1000 384.79678188



THERALASE TECHNOLOGIES INC
LAYER 5 - PWR
PROD: TLC2000 CONTROLLER REV4 DATE: Feb 19 2015 No. of PCB LAYERS: 6



CAM350 V 10.1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: GOV & TO . 1.0 : Eulestions Contact 4DA/EDRA/OCE/DID at CORNER OS PATUS: EULESCOV & TO . 1.0 : EULESC

(b)(4)Proprietary Information - Design	

#### THERALASE TECHNOLOGIES INC

LAYER 6 - SECONDARY SIDE PROD: TLC2000 CONTROLLER REV4 DATE: Feb 19 2015 No. of PCB LAYERS: 6



CAM350 V 10.1.0 : Thus Fab contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188

(b)(4)Proprietary Information - [	Design	

THERALASE TECHNOLOGIES INC	
----------------------------	--

SOLDER MASK - PRIMARY SIDE

PROD: TLC2000 CONTROLLER REV4 DATE: Feb 19 2015 No. of PCB LAYERS: 6

(b)(4)Proprietary Information - Design

THERALASE	TECHNOLOGIES	INC

SOLDER MASK - SECONDARY SIDE

PROD: TLC2000 CONTROLLER REV4 DATE: Feb 19 2015 No. of PCB LAYERS: 6



THERALASE TECHNOLOGIES INC
SOLDER PASTE – PRIMARY SIDE
PROD: TLC2000 CONTROLLER REV4 DATE: Feb 19 2015 No. of PCB LAYERS: 6



CAM350 V 10.1.0 : Thu Feb 19 14:08:52 2015 - (Untitled) : SPP-600-0002-R4.pho Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188



#### THERALASE TECHNOLOGIES INC

SOLDER PASTE - SECONDARY SIDE

PROD: TLC2000 CONTROLLER REV4 DATE: Feb 19 2015 No. of PCB LAYERS: 6

CAM350 V 10.1.0 : Questions? Contact & DA/CORH/QCE/DID-et CORH-FOISTATUS.GOD/@fda/hhs.gov/or-301-796-8188

#### o)(4)Proprietary Information - Design

THERALASE TECHNOLOGIES INC
SILKSCREEN – PRIMARY SIDE

PROD: TLC2000 CONTROLLER REV4 DATE: Feb 19 2015 No. of PCB LAYERS: 6

(b)(4)Proprietary Information -	Design	

	THERALASE	TECHNOLOGIES	INC
- 1			

SILKSCREEN - SECONDARY SIDE

PROD: TLC2000 CONTROLLER REV4 DATE: Feb 19 2015 No. of PCB LAYERS: 6



CAM350 V 10.1.0 : Thu Feb 19 14:08:52 2015 - (Untitled) : NCD-600-0002-R4.drl Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188 (b)(4)Proprietary Information - Design

b)(4)Proprietary Information - Design

# 2.3.6\_Display BlueTooth Schematic

(b)(4)Proprietary Information	- Dəsign						A
							в
							с
(b)(4)Proprietary Informat	ion - Dəsign		Update: 26/ Update: Update:	10/2014	DRAWN CHECKED ENG APPR. Q.A	NAME DATE JACK 26/10/20	<ul> <li>Schematic, TLC2000</li> <li>TFT Display_Bluetooth</li> </ul>
1	2	3	4		5		B 500-0014 R2 SHEET * OF * 6




	SIZE	QTY	SYM	PLATED	TOL
	12	46	+	YES	+/-0.003
	35.43	2	♦	YES	+/003
	39.37	3	×	N0	+/003
a	2				
$\mathcal{D}\mathcal{T}$	ナ				
Ч					

THERALASE TECHNOLOGIES INC
DRILL GRAPHICS
PROD: TLC2000 TFT-Display Bluetooth REV3

DATE: Feb 4, 2015 No. of PCB LAYERS: 2


Refer to : PCB-600-0014-R3.DOC for Fabrication Notes

2.3.7\_Display BlueTooth PCB



THERALASE TECHNOLOGIES INC
LAYER 1 - PRIMARY SIDE
PROD: TLC2000 TFT-Display Bluetooth REV3 DATE: Feb 4, 2015 No. of PCB LAYERS: 2







THERALASE	TECHNOLOGIES INC
LAYER 2	- SECONDARY SIDE

PROD: TLC2000 TFT-Display Bluetooth REV3 DATE: Feb 4, 2015 No. of PCB LAYERS: 2






THERALASE TECHNOLOGIES INC
SOLDER MASK – PRIMARY SIDE
PROD: TLC2000 TFT-Display Bluetooth REV3 DATE: Feb 4, 2015 No. of PCB LAYERS: 2



THERALASE TECHNOLOGIES INC
SOLDER MASK - SECONDARY SIDE
PROD: TLC2000 TFT-Display Bluetooth REV3 DATE: Feb 4, 2015
No. of PCB LAYERS: 2






THERALASE TECHNOLOGIES INC		
SOLDER PASTE - PRIMARY SIDE		
PROD: TLC2000 TFT-Display Bluetooth REV3 DATE: Feb 4, 2015 No. of PCB LAYERS: 2		

CAM350 V 10.1.0 : Wed Feb 04 11:00:59 2015 - (	(Untitled) : SPP-600-0014-R3.pho
--	----------------------------------

(b)(4)Proprietary Information - Design		

THERALASE TECHNOLOGIES INC SOLDER PASTE - SECONDARY SIDE

PROD: TLC2000 TFT-Display Bluetooth REV3 DATE: Feb 4, 2015 No. of PCB LAYERS: 2




SILKSCREEN - PRIMARY SIDE PROD: TLC2000 TFT-Display Bluetooth REV3 PATE: Fab. 4, 2015

 $\rightarrow$ 

DATE: Feb 4, 2015 No. of PCB LAYERS: 2

2.3.7\_Display BlueTooth PCB



THERALASE TECHNOLOGIES INC
SILKSCREEN - SECONDARY SIDE

PROD: TLC2000 TFT-Display Bluetooth REV3 DATE: Feb 4, 2015 No. of PCB LAYERS: 2



2.3.7\_Display BlueTooth PCB



CAM350 V 10.1.0 : Wed Feb 04 11:00:59 2015 - (Untitled) : NCD-600-0014-R3.drl Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188 2.3.8 - Display blue-tooth top picture

(b)(4)Proprietary Information - Design	

# 2.3.9 - Display blue-tooth bottom picture

(b)(4)Proprietary Information - Design	




F			1		
1	2	3	4	5	6
(b)(4)Proprietary Information ·	- Design				
to A the rophotally information	Doolgin				

	2	3	4	5	6
(b)(4)Proprietary Informatio	on - Design				
A					
В					
с					
				NAMES 1 22	
					-2014 TITLE
D				CHECKED	Schematic Laser Diode Driver
				ENG APPR. Q.A	
				Cannot open file C:\Users\Public\Documents\Altium\Al 4\Templates\logo_theralase.hmp	SIZE         NUMBER         REV           D1         B         500-0003-R3         R3
				4\Templates\logo_theralase.bmp	SHEET 2 OF 3
1	2	3	4	5	6


R       1000000000000000000000000000000000000	(b)(4)Pr	oprietary Informatio	² on - Design	۹.	4		t.		6	А
D      CHECKED      CHECKED										
D      CHECKED      CHECKED	в									в
D										
D      CHECKED      CHECKED	с									с
D ENG APPR AND CONTRACT AND CON						-		11-11-2014	TITLE	
1         2         3         4         5         6	D		-				ENG APPR. Q.A Cannot open file C:\Users\Public\Documents\Alt 4\Templates\logo_theralase.bmj	um\AD1	SIZE NUMBER REV B 500-0003-R3 R3 SHEET 3 OF 3	•


### b)(4)Proprietary Information - Design







(b)(4)Proprietary Inform	mation - Design		









7



### (b)(4)Proprietary Information - Design





T



#### (b)(4)Proprietary Information - Design

SOLDER	WASK -	PRINARY	SIDE
PROD: 1	LC2030 90	5/660	DRIVER REV.
DATE: P	OV 21, 20 PCB LAYER	14	





#### b)(4)Proprietary Information - Design







THERALASE TECHNOLOGIES INC SOLDER PASTE - PRIMARY SIDE PROD: TLC2000 905/680 DRIVER REV3 DATE: NOV 21, 2014 No. of PCB LAYERS: 4






#### b)(4)Proprietary Information - Design








#### (b)(4)Proprietary Information - Design








(b)(4)Proprietary Information - Design

## 2.3.13\_Laser Driver Bottom Picture


(b)(4)Proprietary Information - Design

(b)(4)Proprietary Information - Design




 $\bigoplus$ 

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017




 $\odot$ 




## 2.3.15\_Power Pack Pcb

 $\odot$ 






THERALASE TECHNOLOGIES INC	
LAYER 1 - PRIMARY SIDE	
PROD: TLC2000 BATTERY PACK DATE: Nor 10, 2015 No. of PCB LAYERS: 2	(600-0005-83)





THERALASE TECHNOLOGIES	INC
LAYER 2 - SECONDARY	SIDE
PROD: TLC2000 BATTERY DATE: Nor 10, 2015 No. of PCB LAYERS: 2	PACK (600-0005-R3)







THERALASE TECHNOLOGIES INC			
SOLDER MASK - PRIMARY SIDE			
PROD: TLC2000 BATTERY PACK (600-0005-R3) DATE: Nor 10, 2015 No. ef PCB LAYERS: 2			





THERA	ASE TE	CHNOLOGIE	S INC	
SOLE	ER MAS	K – SECON	DARY	SIDE
DATE:	TLC200 Nor 10, PC8 L	2015	PACK	(600-0005-R3)





THERALASE TECHNOLOGIES INC		
SOLDER PASTE - PRIMARY SIDE		
PROD: TLC2000 BATTERY PACK (600-0005-R3) DATE: Nor 10, 2015 No. ef PCB LAYERS: 2		






SILKS	CREEN -	PRIMARY	SIDE	
PROD: TLC2000 BATTERY PACK (600-0005-R3)				
	PC8 LA			

 $\oplus$ 





THERALASE TECHNOLOGIES INC		
SILKSCREEN - SECONDARY SIDE		
PROD: TLC2000 BATTERY PACK (600-0005-R3) DATE: Nor 10, 2015 No. of PCB LAYERS: 2		

(b)(4)Proprietary Information - Design

b)(4)Proprietary Information - Design

# (b)(4)Proprietary Information - Design

### 2.4 Sterilization:

This section of the FDA submittal is not applicable to the Theralase TLC 2000 product as the product is not sold in a sterile condition and is not meant for a one time use. The probe is meant for repeated uses and can be cleaned and disinfected using a cloth moistened with 70% alcohol after each use.

# 2.5 Biocompatibility and Toxicity


(b)(4)Proprietary Information		

# $2.6\ Proposed\ Labeling {\tt ecords}\ {\tt Processed}\ {\tt under}\ {\tt FOIA}\ {\tt Request}\ {\tt 2016-6074}. \ {\tt Released}\ {\tt by}\ {\tt CDRH}\ {\tt on}\ {\tt 1-12-2017}$

b)(4)Proprietary Information

### **Probe Top Label**





Lab Warning Sign TLC-2000

S200-107

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017





# PAIN MANAGEMENT PROTOCOL MANUAL

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188



## **Table of Contents**

CL	INICAL INFORMATION	. 2
	Practical Advice for the Use of Low Level Laser Therapy	. 2
	Treatment Basics	
	Caution	. 2
	General Protocols	. 2
	Acute / Sub-acute / Chronic cases:	. 3
	Indications	3
	Contraindications	3
	Cautions in the use of Laser Therapy	4

### **PROTOCOL INFORMATION**


		_
001 Knee I	Pain	5
oor mice i		-



### **Clinical Information**

### Practical Advice for the Use of Low Level Laser Therapy

Cleaning the skin with isopropyl alcohol will diminish any reflection of the laser beam by the tissue surface due to natural oils.



### **Treatment Basics**

### **General Protocols**


(b)(4)Proprietary Information		

thera

Note:

- Certain pharmaceuticals are light sensitive and laser irradiation may increase or decrease their effectiveness in the body.
- Ask your patients about the sensations they experience during and after treatments. It will provide you with valuable knowledge of treatment protocols and treatment efficacy to apply in future treatments.

### Acute / Sub-acute / Chronic Cases:



### Indications

- 1) All tissue structures: epidermis, dermis, subcutaneous fat, muscle, tendon, bone and ligaments in the knee joint may be treated with a therapeutic laser for bio-stimulative rehabilitative purposes.
- 2) Metal pins, metal plates, plastics and pacemakers are not contraindicated.
- 3) Use clear plastic wrap, cellophane or translucent / transparent bandage types to protect a wound site resident in the knee area prior to laser treatment.

### Contraindications

- 1) Do not treat the abdomen of a pregnant woman.
- 2) Do not treat any types of neoplasia, cancer or cancerous lesions, (benign or malignant).

Copyright © 2015, *Theralase™ Inc.* 




### Cautions in the use of Laser Therapy

- 1) Never point the laser in the direction of any person's eyes.
- 2) Do not treat endocrine glands such as the pancreas and the thyroid.
- 3) For a new patient or a patient who may be light sensitive, spot test the laser on an area of tissue to determine if they are adversely sensitive to laser treatment.
- 4) If a patient is currently prescribed photosensitive medicine, NSAIDS, steroids or corticosteroid medicine, allow 2 weeks as a wash out period prior to laser treatment. If a wash out period is not possible, reduce the laser power to ½ of the recommended power indicated and increase with each treatment as the patient reliance on photosensitive medicine, NSAIDS, steroids or corticosteroid medicine reduces.
- 5) Start treatment at the most painful area first and follow the treatment specifics. Treat surrounding tissue next.





# **Treatment Specification: 001 KNEE PAIN Knee Pain**



### **Treatment Protocol:**

TLC-000 Treatment Record	001 KNEE PAIN
Number of Treatment Points	7

### **Treatment Settings:**

Power	40
Time	1:00
Energy (Joules)	2.4 J

### **Treatment Details:**

b)(4)Proprietary Information

Copyright © 2015, Theralase™ Inc.

### **Treatment Program:**

Frequency	3 per week
<b>Total Visits</b>	12

### 2.8 Advertising and Promotional Material

There are no promotional or advertisement materials available at this time.

### 3.1 Substantial Equivalence

The new product, the TLC-2000 is a Class 3B therapeutic medical laser. It is a new design based on the design of the predicate product, the TLC-1000, also a Class 3B laser. The TLC-1000 product is manufactured and for the past ten years has been sold in Canada under the Canadian Medical Devices License 65816 and has been sold in the United States under a 510(K) Premarket Submission Number K050342.

The TLC-2000 is a more modern approach where the practitioner saves patient diagnosis and treatment protocol information on a tablet computer. (b)(4)Proprietary Information





Predicate Device: TlC-1000

As shown above clockwise from the left:

- 1. The Controller
- 2. The Power Supply
- 3. The Medical AC Cable
- 4. The Probe

New Device: TLC-2000

As shown above clockwise from top left:

- 1. The Tablet Computer
- 2. The power Pack
- 3. The Power Supply
- 4. The Probe
- 5. The Medical AC Cable

(b)(4)Proprietary Information	

### Substantial Equivalence Table

(b)(4)Proprietary Information	

### Substantial Equivalence Table ... Cont'd

(b)(4)Proprietary Information	

The main differences between the predicate Theralase TLC-1000 and the new device, Theralase TLC-2000 Medical Therapeutic Laser Systems are:

1. the additional safeties resulting from the newer laser standard:




2. higher power:

(b)(4)Proprietary Information		

3. the tablet software and database:

(b)(4)Proprietary Information		

4. Power measurement:





Healing at the Speed of Light\* T



Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188

(b)(4)Proprietary Information-Verification Report

April 17, 2015

(b)(4)Proprietary Information-Verification Report

Page 4 of 17

### 2 Software Requirements

o)(4)Proprietary Information-Verification Report

April 17, 2015

Software Validation Report (1).docx





April 17, 2015

Software Validation Report (1).docx

Page 6 of 17

3 Software Verification (b)(4)Proprietary Information-Verification Report	





(b)(4)Proprietary Information-Verification Report	Software Verification and Validation Report	Page 9 of 17
	(b)(4)Proprietary Information-Verification Report	

April 17, 2015

Software Validation Report (1).docx



April 17, 2015

Page 11 of 17

Soltmare + entiteation and + anaat	0
(b)(4)Proprietary Information-Verification Report	
(b)(4)) Tophotary mornation voluidation roport	

Software Verification and Validation Report (b)(4)Proprietary Information-Verification Report	Page 12 of 17
(b)(4)Proprietary Information-Verification Report	

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188

Software Verification and Validation Report	Page 13 of 17
(b)(4)Proprietary Information-Verification Report	

Page 14 of 17

b)(4)Proprietary Information-Verification Report	

April 17, 2015

Page 15 of 17

(b)(4)Proprietary Information-Verification Report	

April 17, 2015

Software Validation Report (1).docx

Page 16 of 17

#### 5 System acceptance test

o)(4)Proprietary Information-Verification Report

April 17, 2015

Software Validation Report (1).docx

1



April 17, 2015

Software Validation Report (1).docx

#### **Risk Management Plan**

Product Description:	Therapeutic Laser
Product Series:	TLC-2000
Date of Analysis:	January 13, 2015
Document Number	<mark>2914</mark>

#### 1) Risk Management Plan

This risk management report is performed on the Theralase TLC-2000 Therapeutic Laser (TLC-2000) for the purpose of investigating the risk of the device. This analysis was planned at the beginning stage of the device design, initial manufacture of device and throughout the life of the device on an annual basis. Main methods and procedures we selected herewith are referred to EN ISO 14971:2007 and EN 60601-1 Edition 3 in regarding the European Medical Device Directives (MDD) on the safety and Performance of the product.

Roger Dumoulin-White, President & CEO and Qing Zhang, Senior Design Engineer have collaborated to compile the risk management plan. Mr. Dumoulin-White is the founder of Theralase Inc. and is well versed in the design, manufacture and application of therapeutic medical lasers. Mr. Zhang is an experienced electronics design engineer who is well versed in the design and manufacture of Theralase therapeutic medical laser systems.

The Theralase TLC-2000 therapeutic medical laser system is designed and manufactured by Theralase Inc. utilizing this risk management plan to minimize risk, maximize safety and efficacy for our customers.

The phases of the life cycle of the Theralase TLC-2000 therapeutic medical laser system covered by the risk management plan include: design, manufacture and post-production usage.

Verification of design changes or process changes to the Theralase TLC-2000 therapeutic medical laser system to minimize risk is completed by Qing Zhang and Roger Dumoulin-White to assess their success in minimizing risk.

The Quality Manager is responsible for all phases of the risk management plan.

Risk that could harm the patient, because bodily injury or death is unacceptable. Risk that prevents operation of the product is acceptable.

#### (b)(4)Proprietary Information-Verification Report

#### 2) Risk Analysis Procedure

a) Intended Use / Intended Purpose and Identification of Characteristics Related to the Safety of the Medical Device

(b)(4)Proprietary Information



3) Estimation of the Risk(s) for each Hazard

See completed FMEA summary below.

#### 3) Risk Evaluation

Risks to be reduced include:

Zero. All risks have been reduced to acceptable levels.

#### 4) Risk Control

#### a) Option Analysis

Risks have been reduced by:

- a. Inherent safety by design
- b. Protective Measures in the medical device itself
- c. Information for safety

#### b) Implementation of Risk Control Measures

All risk control measures have been fully implemented per the attached FMEA summary below.

#### c) Residual Risk Evaluation

No residual risks exist as all risk control measures are in place and all risks have been reduced to acceptable risk levels.

#### d) Risk / Benefit Analysis

Not required. All risks have been reduced to acceptable risk levels.

#### e) Other Generated Hazards

Not required. All risks have been reduced to acceptable risk levels.

#### f) Completeness of Risk Evaluation

Completed per the attached FMEA summary below.

#### 5) Overall Residual Risk Evaluation



FMEA Document TLC D2914 Rev A January 13, 2015

No residual risks exist as all risk control measures are in place and all risks have been reduced to acceptable risk levels.

6) Post Production Information

b)(4)Proprietary Information

7) Failure Mode and Effects Analysis (FMEA) Document Summary

(b)(4)Proprietary Information



Severity Evaluation Criteria:

(b)(4)Proprietary Information

#### 8) Usability Verification Plan and Report

a. Inputs to the USABILITY SPECIFICATION



FMEA Document TLC D2914 Rev A January 13, 2015



**Medical Purpose** 

The TLC-2000 multiple laser probes are used in direct contact with tissue in order to inject photons of light (b) non-invasively and drug-free into tissue. (4)Proprietary



#### Part of body or type of tissue applied to or interacted with

The TLC-2000 multiple laser probes are used in direct contact with tissue depending on what the injury is. Consult protocol manuals for correct placing of the probes on patient.

#### **Intended Operator**

Trained medical practitioner or staff under medical practitioner's guidance and who has taken training on how to use the medical device.

#### Application

Environment:

- Hospital
- Medical Clinic

Conditions of Visibility:

- Ambient Luminance 100 to 500 lux
- Viewing Distance 20 cm to 1 metre
- Viewing Angle orthogonal ± 20°

#### Physical:

- Voltage 100 to 250 vac
- Temperature 0 to 25°C
- 0 to 100 % RH non-condensing

#### (b)(4)Proprietary Information



b. Primary Operating Function



Healing at the Speed of Light\*

b)(4)Proprietary Information

FMEA Document TLC D2914 Rev A January 13, 2015

- Audible indication on Probe to notify user that the laser is on •
- Visible indicator lights to show the power pack is plugged in to the external power supply properly

#### **Resulting Hazardous Situations and Harms**

(b)(4)Proprietary Information		

#### Preliminary review of the User Interface concept

- 1. Interchanging Probes:
  - Probe cable input into probe and power pack clearly identified •
  - Cable connector locks into place and is clearly identified showing location of mating connector •
  - Clear instructions in operator / user manual and quick reference guide Conclusion: No Issues
- 2. Reading Display:
  - Fonts and colours clearly legible from 20 cm to 1 meter ٠
  - Operating instructions for each option clearly identified in Operator / User manual and Protocol manual Conclusion: No Issues
- 3. Setting Up / Preparing for Use:
  - Easy to assemble / disassemble
  - Set-up included in Quick Reference Guide
  - Detailed operating instructions in operator / user manual and on training DVD Conclusion: No Issues
- (b)(4)Proprietary Information 4
- 5. Cleaning:
  - Easy to clean surfaces with Isopropyl Alcohol wipes and swabs ٠
  - Surface finishes are smooth eliminating areas that trap dirt Conclusion: No Issues
- 6. (b)(4)Proprietary Information

d. Use Scenarios

The following use scenario has been performed.



Usability Specification :
Test Description :
Acceptance Criteria :
Test Result :

Internal and / or Customer Specification Reference From Specification clause description of test being undertaken Pass / Fail targets and Pass / Fail criteria Pass or fail is recorded according to the acceptance criteria.

## (b)(4)Proprietary Information



(b)(4)Proprietary Information

Page intentionally left blank.



Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017





## TLC-1000 PAIN MANAGEMENT PROTOCOL MANUAL



## **Table of Contents**

NICAL INFORMATION		 ••••••••••••••••••	
Practical Advice for the Use of Low Le	vel Laser Therapy		
Treatment Basics		 	
Caution			
General Protocols			
Acute Cases:		 	
Sub-Acute Cases:		 	
Chronic Cases:		 	
Treatment Protocol Patient Optimiza	ion	 	
Chronicity:		 	
Darker Skin Colouration:			
Body Mass Index:			
Age:		 	
Fitness level:		 	
Athletes:			
Indications		 	4
Contraindications			

#### **PROTOCOL INFORMATION**

Knee F	Pain	. 6
--------	------	-----

theraLASE

## **Clinical Information**

### Practical Advice for the Use of Low Level Laser Therapy

b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 



b)(4)Proprietary Information

**TLC-1000 THERAPEUTIC MEDICAL LASER** 

Copyright © 2012, *Theralase™ Inc.* 



(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

TLC-2000 Pain Management Protocol Manual Rev 1.0 January 15, 2015

Page 4




(b)(4)Proprietary Information

TLC-1000 THERAPEUTIC MEDICAL LASER



## **Treatment Specification: KNEE PAIN Knee Pain**



#### **Treatment Protocol:**

TLC-2000 Treatment Record	KNEE PAIN
Number of Treatment Points	6

#### (b)(4)Proprietary Information

#### **Treatment Details:**

(b)(4)Proprietary Informatior

#### Copyright © 2012, Theralase™ Inc.

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017





# TLC-1000 OPERATIONS MANUAL



(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



#### TLC-1000 Laser Controller Labels and Symbols

(b)(4)Proprietary Information

July 10, 2012 - Rev. 1.11



(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



## 2. THERALASE TLC-1000 SERIES PRODUCT INFORMATION

b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



(b)(4)Proprietary Information



### The Multiple Laser Probe (TLC-900)

(b)(4)Proprietary Information



(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



(b)(4)Proprietary Information



Hard case for TLC-1000 system (part number: TLC-1200)



# 3. Start-Up Information

<b>b</b>	<b>(4</b> )	Proprieta	rv Infor	mation

Copyright © 2012, *Theralase™ Inc.* 



b)(4)Proprietary Information

unit.

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



Copyright © 2012, *Theralase™ Inc.* 



(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11

thera LA

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017

b)(4)Proprietary Information



b)(4)Proprietary Information

**TLC-1000 THERAPEUTIC MEDICAL LASER** 

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 



(b)(4)Proprietary Information



(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



(b)(4)Proprietary Information



(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



b)(4)Proprietary Information



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 



(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 



(b)(4)Proprietary mormation

Copyright © 2012, *Theralase™ Inc.* 



(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11





Copyright © 2012, *Theralase™ Inc.* 



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 



b)(4)Proprietary Information



## 6. Laser Exposure Level Calculations

b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



(b)(4)Proprietary Information

TLC-1000 THERAPEUTIC MEDICAL LASER

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



(b)(4)Proprietary Information



# 7.Electro Magnetic Compatibility Information

In accordance with clause # 5.2.x of EN 60601-1-2.

Interconnecting Cables

o)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 



(b)(4)Proprietary Information


Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11




(b)(4)Proprietary Information

b)(4)Proprietary Information







Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017





# TLC-2000 OPERATIONS MANUAL

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188



# TABLE OF CONTENTS

# 1. SAFETY INFORMATION – PLEASE REVIEW BEFORE USING THE THERALASE SYSTEM 4

Labelling and an Explanation of Symbols 4 Location of the Labels and Symbols 4 TLC-2000 Laser Controller Labels and Symbols 5 Key Equipment Safety Features 6 Hardware 6 Operational Safety 7 Precautions 8 Safety Warning 8 Disclaimer Information 9

# 2. THERALASE TLC-1000 SERIES PRODUCT INFORMATION 10

Product Overview and Technical Description 10 The Power Supply (TLC-154) 11 The Laser Power Pack (TLC-2000) 12 Laser Output Probes (TLC-2001) 14 Connecting/Disconnecting the Probes to the TLC-2000 Power Pack 15 Accessories 16 Safety Eyewear 16 Patient Safety 18

# 3. START-UP INFORMATION 19

The Power Pack19The Probe20The Tablet22Set-up Procedure23

# 4. SYSTEM OPERATION OVERVIEW 25

Operating the TLC-2000 25 How to Perform a Treatment 25 How to Perform a Measurement 25 Stop Procedure 26 Overview of the Menu 26 System Options Tab 26 Practitioner Tab 27 Administrator Tab 27 System Operation in Detail 27 Common Operations: 27 Turning On/Off the Tablet: 27 Logon to the tablet: 28 Launching the TLC-2000 software: 28 Enter ng Pat ent Informat on 29 Treat ng the Pat ent 30 Ver fy ng the Laser Output 32 Advanced Operat ons 34



How to Create Custom Protoco s 34 How to Check Probe Usage Hours 34 How to Ad ust Probe Buzzer Vo ume 35 Adm n strator Funct ons 35 Add/Remove users 35 Set User Perm ss ons 36 V ew Treatment records 36 Warn ng Messages 37

5

# MAINTENANCE AND SERVICE 38

Inspect on and C ean ng 38 Week y 38 Month y 38 Prevent ve Ma ntenance Schedu e 38 How to check the opt ca output of each aser 38 Transport 38 Calibration/Repair of the Laser Verification Unit 39 Serv ce of Laser D odes 39 Product Warranty 40 Env ronmenta Cons derat ons 40 Thera ase Customer Serv ce Contact Informat on 40

# 6 THERALASE LASER CALCULATIONS

905 nm Ca cu at ons4660 nm Ca cu at ons42Determ ne the Laser C ass43MPE's for sk n and cornea44NOHD Ca cu at on46

# 7 ELECTRO MAGNETIC COMPATIBILITY INFORMATION 49

41



# 1. Safety Information – Please Review Before Using the Theralase System

The Theralase<sup>™</sup> TLC-2000 therapeutic laser is a Class 3B medical laser and as such is a potential hazard for direct and indirect viewing of the laser light.

The Theralase TLC-2000 system as mentioned earlier is composed of four main components; the tablet computer, the Power Supply, the Power Pack and the Laser Probes. The Laser Probes are the applied parts of the system. Safety and safe use of the product are indicated both visually by the use of labels as well as equipment features. These will be detailed below.

NOTE: IT IS IMPORTANT TO BE WELL AWARE OF ALL THESE SAFETY DETAILS PRIOR TO USING THE EQUIPMENT.

NOTE: In Canada, the device is to be installed and operated according to the requirements of the following standard: CSA - Z386-14 Safe Use of Lasers in Health Care.

# Labeling and an Explanation of Symbols





TLC-2000 Laser Probe



TLC-2002 Power Pack

TLC-2003 Power Supp y

F gure Locat on of the safety symbo s on the Thera ase TLC-2000 System



# TLC-2000 LASER CONTROLLER LABELS AND SYMBOLS

b)(4)Proprietary Information



# **KEY EQUIPMENT SAFETY FEATURES**

# HARDWARE

In addition to the safety signs placed on each of the elements of the Theralase<sup>™</sup> system, key safety features are built into the TLC-2000 Power Pack and Tablet software as well. These are illustrated below.



Figure 3. Safety Features

1. Emergency Stop: This button shuts down the entire laser system immediately when pressed.

**2. ON/OFF Key Switch:** This is the main on / off switch for the laser system. When not in use, the keys can be placed in a safe location to prevent unauthorized use of the equipment.

**3. Remote Interlock Connector**: The Remote Interlock Connector is a connector that may be removed for disconnecting the mains power.

**4. Laser Probe:** The Laser Probes are the applied parts of the system. The Probe cannot operate the laser diodes unless a laser protocol has been received from the tablet software.

**5. Tablet Software:** Even if the Key Switch is left in the system with the Remote Interlock Connector in place and the Laser Probe(s) turned on, there is yet another security feature in place. The Laser Probe(s) cannot emit laser light unless treatment data has been received from the Theralase software on the tablet. To login into the operating system the Theralase TLC-2000 application software requires a username and password to operate.





# **OPERATIONAL SAFETY**

Caution: Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure.

b)(4)Proprietary Information

Please note that Theralase<sup>™</sup> supplies Laser Safety Eyewear with all of its products and can be contacted directly at 1-866-The-Lase (1-866-843-5273) for additional sets if required. Order (Part #: TLC-901) directly from Theralase<sup>™</sup> Inc. for all additional laser eyewear needed. The laser eyewear is for diffuse viewing of laser light only.

The Theralase<sup>™</sup> TLC-2000 laser output is athermic. Laser diode package cases may warm up slightly during treatment giving a feeling of warmth to a laser treatment. Under no circumstances should the patient experience a very hot or burning sensation during treatment. If this occurs, discontinue the treatment immediately and contact Theralase<sup>™</sup> Inc. for service.

Observe the safeguards as put forward below. The Theralase™ TLC-2000 therapeutic medical laser is a sophisticated medical device and as such should only be handled by experienced medical practitioners suitably trained in its use.



# PRECAUTIONS

- Make sure that the area to be treated is free from any substance that reflects or prevents penetration of the laser beam. It is recommended that areas to be laser treated are cleansed, prior to treatment.
- Operate the laser in a well-illuminated area, as the pupils of the eye are smaller and not as easily damaged.
- Operate the laser in an area free from mirrors or polished metal surfaces to prevent accidental reflections.

• Operate the laser in a segregated area with a sign posted: b)(4)Proprietary Information

Operate the laser probes in direct contact and perpendicular to the tissue being treated.

Align laser probes to the area to be treated before commencing laser treatment.

# SAFETY WARNING

- The TLC-2002 Laser Power Pack is designed to operate with TLC-2001 Laser Probes only. A maximum of four TLC-2001 Laser Probes may be connected to the TLC-2002 Power Pack at any one time.
- The Theralase™ TLC-2000 laser and all attachments are not designed to be protected from water intrusion and should never be used wet or in the accompaniment of water. If the laser system becomes wet or immersed in water, thoroughly dry the laser system and allow 48 to 72 hours in a warm environment to fully dry prior to use. If problems persist, contact Theralase™ Inc. for service.
- The Theralase<sup>™</sup> TLC-2000 laser and all its attachments should not be used in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide.
- The Theralase<sup>™</sup> TLC-2000 laser and all attachments are not suitable for use in the presence of strong electromagnetic field or other electrical interference.
- The Theralase<sup>™</sup> TLC-2000 laser and all attachments are designed as an integral system and use of parts other than the recognized accessories or parts can degrade the performance of the Theralase<sup>™</sup> TLC-2000 laser and may result in hazardous radiation exposure.



WARNING: No Modification of This Equipment is Allowed WARNING: To Avoid Risk of Electrical Shock, This Equipment Must Only Be Connected to a Supply Mains With Protective Earth.

# DISCLAIMER INFORMATION

The TLC-2000 therapeutic laser system is a medical device and is sold to healthcare practitioners as a tool for the rehabilitation of their patients. The TLC-2000 therapeutic laser system and attachments are utilized to assist in the alleviation of their patient's symptoms. The TLC-2000 therapeutic laser system is not promoted as a cure for any patient ailments.

The information contained in this manual is for reference only and does not preclude the practitioner from using proper judgement in the rehabilitation of their patients.

It is the practitioner's sole responsibility to decide if laser therapy is indicated for their specific patient's rehabilitation.

The treatment protocols detailed in this manual are for reference only and are not promoted as the only protocol available to the practitioner.



# 2. THERALASE TLC-2000 SERIES PRODUCT INFORMATION

The TLC-2000 Therapeutic Laser System is used in the adjunctive treatment of knee pain.

The TLC-2000 laser probes are used in direct contact with tissue in order to inject photons of light (b)(4)Proprietary (b)(4)Proprietary

### PRODUCT OVERVIEW AND TECHNICAL DESCRIPTION

The Theralase™ TLC-2000 series of products is comprised of four main components

- The Tablet Computer with Theralase software (p/n TLC-2004)
- The Power Supply (p/n TLC-2003)
- The Laser Probe(s) (p/n TLC-2001)
- The Power Pack (p/n TLC-2002)

The Power Pack only supplies power to the Laser Probe(s). The Laser Probes are controlled by Theralase software resident on the Tablet. Under software control the Laser Probes operate in a patient treatment mode with the power level and time based on the protocol resident in the tablet software. Communications between the Tablet and the Laser Probe(s) is via a wireless Blue Tooth® link (b)(4) Proprietary Information

b)(4)Proprietary Information

# (b)(4)Proprietary Information



(b)(4)Proprietary Information



# THE LASER POWER PACK (TLC-2002)

The Theralase<sup>™</sup> TLC-2000 laser system Power Pack (TLC-2002) is shown below. The Power Pack consists of an aluminum case with a plastic top and bottom. There is a separate plastic stand that the bottom of the Power Pack fits into. There are four ports to accommodate up to four TLC-2001 Laser Probes.



Laser Probe Port



Unit ON Indicator



Laser Probe Power ON Indicator



Remote Interlock Connector



Emergency (E-STOP) Switch



Key Switch

Figure 5. The TLC-2000 system Power Pack (TLC-2002)

The Theralase<sup>™</sup> TLC-2000 laser system and all attachments are rated as Class 1 ME equipment, type BF and are equipped with Class 3B lasers.

Details regarding the safe use of the TLC-2000 system components are detailed in section 2 and operation of the unit in sections 3 and 4.

The top of the unit has an emergency STOP button (red button with white arrows), a key switch for switching the unit off and being able to lock the unit from unauthorized use, an interlock switch and a green ON LED. A further detail of these items as well as their operation is detailed in section 2.

The following are the detachable parts for the TLC-2000:

Power Supply:TLC-2003 Power Supply (+24 VDC, 4.2 A)Laser Probe:TLC-2001 Laser Probe



# WARNING: PLEASE DISCONNECT POWER SUPPLY BEFORE SERVICING

The Power Supply (TLC-2003)

The TLC-2003 is a universal power supply – the input voltage range is 100VAC to 240VAC, 47 to 63 Hz, meaning that it can be safely used anywhere in the world. The TLC-2003 power supply is able to power the TLC-2002 Power Pack.



Plugs in to main AC line via power cord (TLC-155)

Power output: Voltage: +24 VDC Current: 4.0 A



Power supply TLC-2003 plugs in to the Power Pack - TLC-2002

Figure 7. TLC-2003 Power Supply and associated Line Cord

With the proper cord set the TLC-2003 can be plugged in anywhere in the world. The cord set (TLC-155) supplied with the Theralase™ TLC-2003 power supply is a hospital grade 115 to 120 VAC, 50 to 60 Hz approved cord set. This is built to a North American standard and was chosen to work in Canada and the United States.

Some of the key features of the TLC-2003 Power Supply include

100-240 VAC Universal Input In Case IEC320 Single Output up to 101 Watts Outputs Regulated with Low Ripple Meets Safety Agency Requirements Complies with EMC/EMI Regulations CE Compliant Impact Resistant Polycarbonate Enclosure

# CAUTION:

• Use only the TLC-2003 Power Supply (+24 VDC, 4.2 A) to power the TLC-2002 Power Pack. This is the only Power Supply approved to be used with the TLC-2000 laser system.

• Outside of North America, the proper international cord set must be specified and purchased from Theralase<sup>®</sup> Inc. or purchased separately for the country of preference.



LASER PROBES (TLC-2001) b)(4)Proprietary Information

Figure 10. Top and Bottom view of the Laser Probe (TLC-2001)

(b)(4)Proprietary Information

CAUTION

The Laser Probe Cover should be placed over the lens area when not in use.



# Connecting/Disconnecting the Laser Probes (TLC-2001) to the Power Pack (TLC-2002)

All the Laser Probes feature a very simple medical grade push-pull cable connector. This consists of a cable with a plug assembly on each end and a printed circuit board-mount receptacle assembly in the Power Pack and the Laser Probe.

Assurance of proper mating is provided by keying. The key(s) of the cable plug assembly are designated by an arrow on the plug and must align with the key of the receptacle assembly.




Figure 11. Probe Assembly illustrating the key

The cable plug assembly must be properly aligned with the receptacle assembly, and then inserted until the latches engage the receptacle assembly. There should be no need to force the connector into the receptacle. To ensure proper mating, there must be an audible "click". Refer to Figure 12.



Figure 12. Mating the Connectors - Cable Plug Assembly Key is aligned with Receptacle key

The connectors must be unmated by first gently sliding the housing subassembly away from the receptacle assembly – see Figure 13 and then pulling the entire connector out. PLEASE hold the main body and DO NOT pull the cable.




Figure 13. Un-mating the connected assembly



# ACCESSORIES

Accompanying the TLC-2000 system, in addition to the above main items, will be two sets of laser safety goggles (TLC-901).

# Safety Eyewear:

Two sets of laser safety goggles (part #: TLC-901) are supplied by Theralase as standard accessories. These laser safety goggles are provided for viewing of diffused laser light only and are not indicated for direct viewing of the laser output.

It is highly recommended that both the practitioner and the patient wear laser safety eyewear during laser therapy treatment.



Figure 14. Theralase supplied Laser Safety Eyewear

Optical Density ("O.D.") is the scientific term to describe a material's ability to block light of a certain wavelength. It is a logarithmic formula, where the higher the O.D. value, the better the material blocks light. For example, an O.D. of 1 means the material will allow 10% of the light penetrate the eyewear at a particular wavelength.

	O.D.	Light Penetration	Light Absorption
(b)(4	l)Proprieta	ry Information	



(b)(4)Proprietary Information

# CAUTION:

Use only Theralase<sup>™</sup> Inc. supplied laser safety eyewear or laser safety eyewear that meets or surpasses the above noted optical densities at the above-specified wavelengths

NEVER POINT A LASER TOWARDS A PERSON'S EYES UNDER ANY CIRCUMSTANCES



Figure 14. Hard case for TLC-2000 system (part number: TLC-1200)



# CONTRAINDICATIONS IN THE USE OF LASER THERAPY:

- 1. Do not treat the addomen of a pregnant woman.
- 2. Do not treat any tpes of neoplasia, cancer or cancerous lesions (benign or malignant).

# CAUTIONS IN THE USE OF LASER THERAPY:

- 1. Never point the laser in the direction of any person's eyes.
- 2. Do not treat endocrine glands such as the pancreas and the thyroid.
- 3. For a new patient or a patient who may be light sensitive, spot test the laser on an area of tissue to determine if they are adversely sensitive to laser treatment.
- 4. If a patient is currently prescribed photosensitive medicine, NSAIDS, steroids or corticosteroid medicine, try to allow 2 weeks as a wash out period prior to laser treatment. If a wash out period is not possible, reduce the laser power to 1/2 of the recommended power indicated and increase with each treatment as the patient reliance on photosensitive, NSAIDS, steroids or corticosteroid medicine reduces.
- 5. Start treatment at the most painful area first and follow the treatment specifics. Treat surrounding tissue next.
- 6. Under no circumstances should the patient feel anything more than a slight warming effect. If the patient feels hot sensations, pain or burning, STOP LASER TREATMENT IMMEDIATELY and check laser equipment. If laser probe head is hot to the touch, call immediately for service. Detailed contact information is provided in Section 5 Maintenance and Service, of this manual.
- 7. Under no circumstances over-treat an individual as over-treatment may lead to bio-inhibition of tissue structures.



# 3. START-UP INFORMATION

A detailed description of the individual system components

# The Power Pack (TLC-2002)

The TLC-2002 Power Pack provides power to each TLC-2001 Laser Probe. The probes ports, numbered 1 to 4 are the same electrical output. A TLC-2001 Laser probe should always be connected to location #1 first, then #2, etcetera.



Figure 16. A detailed look at the Power Pack

**1. Remote Interlock Connector** The Remote Interlock Connector is capable of disconnecting the TLC-2001 Laser Probe(s) from main power when removed. To remove, press the button on the connector and pull the Remote Interlock Connector straight out. Main power is now disabled and no action will restore it until the Remote Interlock Connector is replaced. To replace, insert straight into the matching receptacle.

2. Emergency Stop This button shuts down the entire laser system immediately when pressed, disabling mains power.

# NOTE: Depressing the Red Emergency Stop Button immediately disables the laser system. Operation may be resumed by rotating the Red Emergency Stop Button 1/4 of a turn to the right and allowing it to "pop out" again.

**3. Key Switch** This is the main on / off switch for the laser system. When the key is turned to the left, the system is off, disabling mains power. When the key is turned ¼ turn to the right, the system is on. The key can only be removed when in the OFF position.



TLC-2003 Power Supply input



Connection to the TLC-2001



**4. Power Supply Input** Input for black DC adapter cord from the TLC-2003 Power Supply. When the DC adapter is connected the TLC-2000 Power Pack may be operated directly from an AC wall outlet.

Caution: Use only the TLC-2003 Power Supply for the Theralase<sup>™</sup> system.

**5. Probe Power Ports** Power supply ports for up to four (4) TLC-2001 Laser Probes (High Density 6 pin medical grade connector).

**6. Power IN LED** This LED is illuminated only when the external power supply is connected and the key switch is in the ON position, the Remote Interlock Connector is installed and the Emergency Stop is not activated.

7. Power OUT LED This LED is illuminated only when both of the internal power supply voltages to the probes are good.

Note: The TLC-2001 Laser Probes will not function if this LED is OFF.

# The PROBE

b)(4)Proprietary Information



# b)(4)Proprietary Information

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188



# THE TABLET (TLC-2004)

The Tablet comes with the Theralase® software already installed and the keyboard and touch screen together. The touch screen can be separated from the keyboard if desired. There is a pop up touch screen keyboard; however, it may be easier to input data using the attached keyboard. With the keyboard attached, selections can be made using the cursor or the touch screen.

When delivered, the Tablet has an administrator username and password. It is recommended that the system be configured for new passwords and user names. Refer to the section Setting Username and Passwords in the Administrator Menu section for this procedure.

The Tablet is shipped with the manufacturer's user manual and software license card. Refer to those documents for additional information.



Figure 16. A detailed look at the tablet

- **1. Power On** Hold this button down for two (2) seconds until the touch screen turns on.
- **2. Charging Port** Plug the micro USB connector from the charger into this port to power the tablet from AC voltage.
- 3. USB Port Plug the Blue Tooth<sup>®</sup> USB stick into this port.
- **4. Release Button** Press this button to separate the keyboard and the touch screen.

**5. Touch Pad** Use this pad to control the cursor and make selections by moving your finger over the surface of the touch pad. Press the lower left corner of the touch pad to make a selection.



# SET-UP PROCEDURE

1. Unpack the contents of the case and identify each of the pieces as outlined in the figure below.

Identify the following components of your TLC-2000 system:

- 1. Tablet computer (TLC-2004)
- 2. Tablet power supply and cable
- 3. Power Supply (TLC-2003)
- 4. Medical AC power cable
- 5. Power Pack (TLC-2002)
- 6. Laser Probe(s) 1 to 4 (TLC-2001)
- 7. Probe cable(s) 1 to 4
  - 8. Straps
- 9. Safety Eyewear



Figure 20. TLC-2000 system components

Set the TLC-2002 Power Pack on a steady surface and ensure the key switch is in the "off" position. Make all connections as shown below and read the following set of instructions BEFORE operating the system.



Figure 20. How the TLC-2000 system connects together



2. Connect the 120 VAC hospital grade line cord (TLC-155) into the TLC-2003 Power Supply (or approved Hospital Grade cord set for country of use) and then connect the TLC-2003 power supply into the power input located on the Power Pack.

3. Connect the Laser Probe(s) (TLC-2001) to the Power Pack (TLC-2002) by inserting the medical grade push-pull connector to mate with the receptacle in the Laser Probe port.

**Important Note:** Port 1 **MUST** be used for one of the Laser Probes. Detailed instructions on how to connect the Laser Probe to the Power pack is given in section 3. Leave the Laser Probes turned OFF.

Note: Connect the Laser Probe(s) to the Power Pack PRIOR to switching the Power Pack ON.

4. Connect the Tablet power supply to the mirco USB connector on the Tablet and plug into a 120 VAC source. The power supply is universal and will need a plug adaptor for 240 VAC.

5. Open the Tablet and turn it on. Then log on to the operating system with the supplied username and password.

6. Select the TLC-2000 icon and enter the username and password to start the application.

Note: There is a user name and password for the Tablet computer and also for the TLC-2000 software.

- 7. At this point the Welcome screen should be displayed with icons for the actions.
- 8. Turn on the Power Pack by:
  - a. Turning the Emergency Stop switch a quarter turn clockwise to allow it to "pop" out
  - b. Inserting the Remote Interlock Connector into its receptacle on the front panel of the Power Pack until it locks firmly into place
  - c. Turn the Key Switch 90° to the right
- 9. At this point the Power Pack Power IN and Power OUT LEDs should illuminate.
- 10. The system components are now set up and ready to be used.

# **REMINDER:**

The Theralase software must be operating with patient data in the database **BEFORE** the Laser Probes will emit laser light.



# **4. SYSTEM OPERATION OVERVIEW**

o)(4)Proprietary Information



# **Stop Procedure**

The Laser Probe may be turned off during operations by either of the following means:

- 1. Slide the switch on the Laser Probe to the OFF position.
- 2. Press the STOP button on the Tablet software
- Press the Red Emergency Stop Button to stop Laser Probe operation and disable the mains power to the TLC-2000 laser system. Laser operation can only be re-started by turning the E-STOP switch a quarter turn clockwise to allow it to "pop" out.
- 4. Turn the on / off power key switch to the left ¼ turn to stop laser operation and to disable the mains power to the TLC-2000 laser system.

The Key Switch or Remote Interlock Connector may now be removed if desired to disable the laser system from unauthorized operation.

**Note:** Even if the Key Switch is left in the system and the Remote Interlock Connector is in place, the security password must be known to enable laser operation; thereby, providing a 3 step security protocol to prohibit laser use by unauthorized personnel.

# **Overview of the Menu**







PRACTITIONER TAB



(b)(4)Proprietary

**ADMINISTRATOR TAB** 

b)(4)Proprietary Information

# System Operation in Detail

COMMON OPERATIONS

# **Turning On the Tablet:**

Turn on the tablet by pressing and holding the power on button on the lid of the tablet for two (2) seconds. Keep it pressed until the screen starts.



### **Turning the Tablet OFF**

To turn the tablet OFF, swipe the screen in a Right to Left motion from the right hand edge of the screen. A vertical menu will be displayed. (To close the vertical menu, swipe the screen in a Left to Right motion.) Touch the Settings Icon to display the settings menu. Touch the power icon and select SHUT DOWN from the pop up window.



(b)(4)Proprietary Information



# b)(4)Proprietary Information

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188



b)(4)Proprietary Information



# (b)(4)Proprietary Information

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188



b)(4)Proprietary Information



(b)(4)Proprietary Information



# **Advanced Operations**

#### **How to Create Custom Protocols**

b)(4)Proprietary Information







# WARNING MESSAGES



# 5. Maintenance and Service

This section deals with maintenance and serviceability of the Theralase<sup>™</sup> system by the end user. If further detailed service of your laser system is required, contact Theralase Inc directly for assistance. Detailed contact information is supplied at the end of this section.

#### **Inspection and Cleaning**

PRE AND POST TREATMENT

WEEKLY

Follow the stop procedure in the basic operation section of this manual. Remove the key from the power pack. Disconnect the Laser Probe(s) and power input to the Power Pack. Carefully clean each of the TLC-2001 Laser Probes with isopropyl alcohol swabs being careful not to immerse the Laser Probe heads. Remove the lense cap and immerse a Q-tip into isopropyl alcohol and rotate in the laser aperture carefully to avoid damaging the laser diodes.

#### MONTHLY

Follow the stop procedure in the basic operation section of this manual. Remove the key from the power pack. Disconnect the Laser Probe(s) and power input to the power pack. Inspect the laser diodes for cracks or chips in the glass. A visible chip in the glass, a crack in the glass or a missing glass mandates replacement of the laser diode. Contact Theralase immediately for replacement of the defective / damaged laser diode.

#### PREVENTIVE MAINTENANCE SCHEDULE



#### Transport

Remove the key from the main power switch and store in a safe location. Disassemble the laser system, carefully packing all components and taking special care with the laser probe assembly. In packing and unpacking the apparatus, be careful not to twist or kink the probe cable(s) too sharply to prevent damage to the probe cable assembly. Ensure proper storage of the Laser Probesto avoid laser diode damage. If you damage your laser diodes, please notify Theralase<sup>™</sup> Inc. immediately for service / replacement. (The laser system should be kept within 0 to 40°C, 500 hPa to 1068 hPa and 0 to 80% relative humidity during transport).



#### Caution:

Calibration must only be performed by Theralase<sup>™</sup> Inc. or by trained personnel possessing the proper knowledge and test equipment. Contact Theralase<sup>™</sup> Inc., if after testing, software indicates maintenance of the laser diodes is required.

#### Calibration/Repair of the Laser Verification Unit

Calibration or repair of the real-time Laser Probe Cover is not a process that can be performed in the field. There are no user serviceable parts in the probe. In the event, if such a situation does arise, it is recommended to call Theralase Customer Support at 1-866-843-5273 or email service@theralase.com for detailed instructions on how to best handle this condition.

#### Service of Laser Diodes

The TLC-2000 therapeutic medical laser system contains no user serviceable parts. If the laser diodes require maintenance then contact the factory and return the system for service.



Product Warranty

(b)(4)Proprietary Information

# **Environmental Considerations**

This product contains electronic components and as such is considered e-waste and must be disposed of responsibly. Do not dispose of this product with waste destined for landfill or incineration. This product must be recycled according to local laws & regulations. At the end of this product service life, estimated to be 10 years, the product can be returned to Theralase for disposal.

# **Theralase Customer Service Contact Information**

Please note the following contact info for Theralase Customer service. When sending your system in, it would be advisable to make a copy of your checklist (see Section 2) and attach this to your package. Describe the issue that your therapeutic laser system is experiencing and the customer service department will instruct you on the recommended service / re-placement parts for your equipment.



























#### 510(k) Safety and Efficacy

Theralase Inc.

TLC-2000 Therapeutic Medical Laser System

June 19, 2015

#### **General Information** I.

<b>Device Generic Name:</b>	TLC-2000 Laser
Trade Name:	TLC-2000 Therapeutic Medical Laser System
Device Classification:	Class II, Performance Standards 21CFR Part 890.5500 Infrared Lamp, Non-heating
Product Code:	NHN
Applicant Name and Address	: Theralase Inc. 1945 Queen Street East Toronto, Ontario M4L 1H7 Canada Phone: (1) 416.699.5273 Fax: (1) 416.699.2250
510(k) Number:	Pending

#### **II.** Device Description:

The TLC-2000 Therapeutic Medical Laser System is a non-heating lamp as defined in 21 CFR 890.5500. The TLC-2000 System is a precision medical device consisting of a hand held portable multiple laser probe (b)(4)Proprietary Information

a tablet computer and a power supply.

## **III. Indications for Use:**

The TLC-2000 Therapeutic Medical Laser System is indicated for adjunctive use in the temporary relief of pain associated with knee disorders.

## **IV. Contraindications and Warnings:**

- A. Contraindications:
  - Direct irradiation of the eyes Class 3B lasers are potentially harmful to the retina, although retinal damage is highly unlikely from exposure to the laser beam of the TLC-2000 Laser when used as directed in the Operator's Manual and under normal use by a licensed healthcare practitioner.
  - Pregnancy The Laser is contraindicated for use in the near vicinity or over the pregnant uterus.
  - The TLC-2000 is a symptomatic treatment, and its use is contraindicated where analgesia may mask progressive pathology, and should not be used where the practitioner would normally avoid the use of any analgesia in order to retain the beneficial aspects of pain.
- B. Warnings
  - Refer to the warning messages in the TLC-2000 Therapeutic Medical Laser Operator's Manual in Appendix C.

## **V. Alternative Procedures or Practices**

The primary alternatives for the use in treating chronic or acute knee pain include the use of TENS, steroids, non-steroidal anti-inflammatory drugs (NSAIDs) and barbiturated drugs (known to be addictive and relatively toxic).

## VI. Adverse Events

There are no known safety or adverse effects reported in the literature from the use of Low Level Laser Therapy. There were no adverse effects reported by the principal investigators during the course of the clinical study that was performed for the predicate device.

The potential incidents, although the likelihood of occurrence is extremely rare, may include: eye pain, retina or other ophthalmic damage if the laser beam is pointed directly at or into the eye for prolonged periods of time; transient redness at the treatment site; swelling at or near the area of laser irradiation; tingling or itching and/or temporary paresthesia.

## **VII. Efficacy**

The TLC-2000 Therapeutic Medical Laser System utilizes the very same protocol, for the temporary relief of pain associated with knee disorders, as the predicate device and therefore the clinical study results for the predicate device can be applied to the new device.

11 \ / 4	<b>D</b> •			
(b)(4)	)Proprie	tarv In	format	lon
		tary in	format	IUII

10 2017

2016

6074 D

D

### Annex A (informative)

# **Biological evaluation tests**

# **Engineering Test Report**

	Engineering rest report	
(	b)(4)Proprietary Information	



(b)(4)Proprietary Information	



(b)(4)Pro	prietary	Intorma	tion
	priotary	innonna	lion

## Wavelength Measurements

(b)(4)Proprietary Information		

Error - Couldn't merge file with following reason - PdfReader not opened with owner password 09002621826d0372.pdf

System attempted to attach the file. Please look at attachments to open this file manually.

HEALING AT THE SPEED OF LIGHT\*

## **TLC-2000 Thermal Tests**


June 24, 2015

)(4)Proprietary Inform	ation-Testing		

From the data file the maximum temperature recorded was 37.8<sup>0</sup> Celsius on the laser diode case. The external case reached a maximum temperature of 36.5<sup>0</sup> Celsius.

Data File:

Header Size: 8							
Version: 2							
Sampling Interval: 2							
Sampling Rate: 0.5							
Sample Count: 1234							
Device Serial Number:							
0							
Culture Info: en-US							
Sample Number	Date/Time	<b>CHANNELO</b>	CHANNEL1	CHANNEL2	CHANNEL3	Events	
	ropriotory	rmation Tastin	a			DAQ	
(D)(4)P	rophetary inic	rmation-Testin	y			art	




Page 4

R	ecords Process	ed under FOIA	Request 2016-60	74. Released by Cl	DRH on 1-12-2017
2)(Z	1)Proprietary	Information	-Testing		

Page 5




.




Page 7

Page 9











Page 20

Page ZI

b)(4)Proprieta	ry Informatior	r-Testing

(b)(A) Duran vistams in formation. Testing
(b)(4)Proprietary Information-Testing

HEALING AT THE SPEED OF LIGHT\*

## TLC-2000 Skin Temperature versus Power Level

b)(4)Proprietary Information-Testing	
and the photon of the stand of	

#### Specification


(b)(4)Proprietary Information-Testing

#### Specification

(b)(4)Proprietary Information-Testing

**Test Results** 

b)(4)Proprietary Information-Testing



Observations: b)(4)Proprietary Information-Testing

#### Specification

Conclusion:

(b)(4)Proprietary Information-Testing

## **Usability Evaluation Report**

## **TLC-2000 Therapeutic Laser System**



Questions? Contact

r 301-796-8188

### **Table of Contents**

#### 1.0 Medical Device

(b)(4)	Proprie	tarv Inf	formatio	n-Test	tina
$(\mathbf{r})(\mathbf{r})$	порне	tary m	ormatio	n- I esi	ung

#### 1.3 Medical Device Description

(b)(4)Proprietary Information-Testing	

#### Laser Treatment Probe (TLC-2001):

(b)(4	)Proprie	tarv Inforn	nation-Te	stina
	л төрнө	aly mom	Ialion-Te	Sung



Tablet (TLC-2004):

Power Pack (TLC-2002):

(b)(4	Prop	riotary			- Loet	
(b)(4	)Prop	rietarv	Inforr	nation	- I OST	na

#### The Power Supply (TLC-2003):

(b)	(4)	P	ro		01		In						0			
D				ori		$\sim$	In	ГО	m	na	Ο	n-	-	51	In	C

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188

Specifications		
Specifications (b)(4)Proprietary Info	rmation-Testing	

#### 1.4 Intended Medical Use

#### 2.0 DUI

(b)(4)Proprietary Information-Testing		
(b)(4)Fropholary monnation-resulty		

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188

#### 4.0 DUI Design

#### 5.0 Evaluation Planning

(b)(4)Proprietary Information-Testing

6.0 Validation Testing b)(4)Proprietary Information-Testing





Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017



U.S. Food and Drug Administration Center for Devices and Radiological Health Document Control Center - WO66-G609 10903 New Hampshire Avenue Silver Spring, MD 20993-0002 FDA CDRH DMC OCT 3 0 2015 Received October 29, 2015 KISI8|6/SOOI

(b)(4)Proprietary Information

1945 Queen Street East, Toronto, Ontario M4L 1H7, Canada •416.699.5273 Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188 Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017



(b)(4)Proprietary Information

6

1945 Queen Street East, Toronto, Ontario M4L 1H7, Canada •416.699.5273 Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@tda.hhs.gov or 301-796-8188 Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017



(b)(4)Proprietary Information

ł.

1945 Queen Street East, Toronto, Ontario M4L 1H7, Canada •416.699.5273 Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@tda.hhs.gov or 301-796-8188

K15181615001

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017

# HEALING AT THE SPEED OF LIGHT®

U.S. Food and Drug Administration Center for Devices and Radiological Health Document Control Center - WO66-G609 10903 New Hampshire Avenue Silver Spring, MD 20993-0002

#### Re Submission: K151816

b)(4)Proprietary Information

October 27, 2015

FDA/CDRH/DCC OCT 2 9 2015 RECEIVED

1945 Queen Street East, Toronto, Ontario M4L 1H7, Canada Ourestions 3 Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 3011796-8188 •410.099.5273 •1.866.845.5273 •WWW.theralase.com Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017



(b)(4)Proprietary Information

1945 Queen Street East, Toronto, Ontario M4L 1H7, Canada Ouestions? Contact FDA/CDRH/OCE/DID at CDRH FOISTATUS.GOV@fda.hhs.gov or 301796-8188 •416.699.5273 •1.866.843.5273 •WWW.theralase.com



Please refer to the attached revised documents:

(b)(4)Proprietary Information

1945 Queen Street East, Toronto, Ontario M4L 1H7, Canada •415.699.5273 •1.866.843.5273 •WWW.theralase.com Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017



October 29, 2015

U.S. Food and Drug Administration Center for Devices and Radiological Health Document Control Center - WO66-G609 10903 New Hampshire Avenue Silver Spring, MD 20993-0002

b)(4)Proprietary Information

1945 Queen Street East, Toronto, Ontario M4L 1H7, Canada 1945 Gon Street East, Toronto, Ontario M4L 1H7, Canada



018\_Instructions for Use 020\_Substantial Equivalence 024\_A\_TLC2000 Operation Manual 0)(4)Proprietary Information

1945 Queen Street East, Ioronto, Ontario M4L 1H7, Canada



#### Please refer to the attached revised documents:

b)(4)Proprietary Information

1945 Queen Street East, Toronto, Ontario M4L 1H7, Canada 1945 Gonga FBA/CDRH/0012/0006 CDRH3F508743 US.GOV@100.0012/0012001/200818801 DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

#### Indications for Use

510(k) Number *(if known)* K151816

Device Name

TLC-2000 Therapeutic Medical Laser System

Indications for Use (Describe)

The Theralase TLC-2000 Therapeutic Medical Laser System is indicated for adjunctive use in the temporary relief of pain associated with knee disorders with standard chiropractic practice.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

#### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov* 

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Form Approved: OMB No. 0910-0120 Expiration Date: January 31, 2017 See PRA Statement below.

#### **1.8 Product Description**

The Theralase<sup>™</sup> TLC-2000 Therapeutic Medical Laser System is a Class 3B medical laser system and is indicated for "Adjunctive Use in the Temporary Relief of Pain Associated with Knee Disorders with Standard Chiropractic Practice." The Theralase<sup>™</sup> TLC-2000 Therapeutic Medical Laser System is comprised of the following components:

(b)(4)Proprietary Information	



The external power supply is a medically approved power supply. The input voltage range is 100 240 VAC with a frequency range of 47 73 Hz. The power supply connects to the AC source with a medically approved ac cord.

The system is supplies with the TLC-2000 Operations Manual, Pain Management Protocol Manual, safety glasses and a carrying case.

### **2.1 Device Description**


#### **System Component Description:**

#### The Probe:



- 1 Display
- 2 Display Window
- 3 Power On/Off Switch
- 4 Start/Stop Button
- 5 Laser Diode Aperture
- 6 Optodes
- 7 Probe Body
- 8 Strap Handles
- 9 Laser Output Window

The probe connects to the power pack with an 8 foot long cable. The probe has a color Liquid Crystal Display, power on/off switch and a treatment start/stop switch. The probe contains the b x x x must have a long the circuitry to operate them. There is an



#### The Tablet:




- 1 Power On
- 2 Charging Port
- 3 USB Port
- 4 Release Button
- 5 Touch Pad

The tablet is a windows based computer running Windows 8.1 operating system. The tablet has a Blue Tooth 4.0 radio and WiFi. The tablet has a 10.1" touch screen and a detachable keyboard. The Theralase application software is loaded on the tablet when supplied with the system.

The Theralase TLC-2000 application software is used by the practitioner to enter patient data, to select the knee pain treatment protocol, which is indicated for Adjunctive Use in the Temporary Relief of Pain Associated with Knee Disorders with Standard Chiropractic Practice, to store patient treatment data and to transmit the treatment settings to the probe. For all data that is input, the software will display a verification window listing the data and ask the operator to confirm that all the data is correct before saving to the database.

#### **The Power Pack:**



- 1 Remote Interlock
- 2 Emergency Stop
- 3 Key Switch
- 4 Power Input
- 5 Probe Power Ports
- 6 Power In / Power Out LED indicators

The power pack contains the electronics which provide the operating voltages for the probe. The power pack contains the necessary safety systems in the form of an emergency stop switch, a remote interlock and a key operated on/off switch.

#### The Power Supply:

The power supply is a medically approved universal, 100 240 VAC input voltage rating at 47/63 Hz, power supply. The output is 24 volts DC for the power pack. A medically approved AC power cord is included.



#### b) ACCESSORY COMPONENT DESCRIPTIONS

#### Safety eyewear:

Supplied with the system are two sets of standard, TLC-901 laser safety glasses one pair for the clinician and one pair for the patient. These laser safety glasses are provided for viewing diffused laser light only and are not for direct viewing of the laser output.





## c) Specifications

(b)(4)Proprietary Information		

#### d) Intended Use of the TLC-2000

The TLC-2000 therapeutic medical laser system is indicated for Adjunctive Use in the Temporary Relief of Pain Associated with Knee Disorders with Standard Chiropractic Practice.

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188 Label Set TLC2000 P/N 701-0021-R1 Material Specification - UL approved xxxxxxx

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017





# TLC-2000 PAIN MANAGEMENT PROTOCOL MANUAL

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188



## **Table of Contents**

Clinical Information	. 2
Practical Advice for the use of Low Level Laser Therapy	2
Cautions in the use of Laser Therapy	. 2
Oper ation Procedure (Knee Protocol)	. 3
Treatment Specification: 001 KNEE PAIN	.4



### **Clinical Information**

Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.




#### **Operation Procedure (Knee Protocol)**

Caution: The patient and the practitioner should always wear protective eyewear during laser operation.



## **Treatment Specification: 001 KNEE PAIN** Adjunctive Use in the Temporary Relief of Pain Associated with Knee Disorders with Standard Chiropractic Practice.



#### Figure 1

#### **Treatment Protocol:**

TLC-000 Treatment Record	001 KNEE PAIN
Number of Treatment Points	7





#### 3.1 Substantial Equivalence

The new product, the TLC 2000 is a Class 3B therapeutic medical laser. It is a new design based on the design of the predicate product, the TLC 1000, also a Class 3B laser. The 1000 product is manufactured and for the past ten years has been sold in Canada under the Canadian Medical De ices License 65816 and has been sold in the United States under a 510(K) Premarket Submission Number K050342.


The TLC 2000 is a more modern approach where the practitioner sales patient treatment protocol information on a tablet computer. (b)(4)Proprietary Information (4)Proprietary Information

(b)(4) Bropriston (lafora	Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017
(b)(4)Proprietary Inforn	

## Substantial Equivalence Table

h)///)Proprietory Information	_
b)(4)Proprietary Information	

## Substantial Equivalence Table ... Cont'd (b)(4)Proprietary Information



Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017





## TLC-1000 PAIN MANAGEMENT PROTOCOL MANUAL

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188



### **Table of Contents**

CLINICAL INF	ORMATION	2
Practical A	Advice for the Use of Low Level Laser Therapy	2
Treatmen	t Basics	2
	on	
Gener	ral Protocols	2
Indication	18	3
Contraind	lications	3
Cautions i	in the use of Laser Therapy	4

#### **PROTOCOL INFORMATION**

001 Knee Pain
---------------

theraLA 

## **Clinical Information**

#### Practical Advice for the Use of Low Level Laser Therapy

b)(4)Proprietary Information


Copyright © 2012, *Theralase™ Inc.* 



#### Cautions in the use of Laser Therapy




## Treatment Specification: 001 KNEE PAIN Knee Pain



#### **Treatment Protocol:**


TLC-000 Treatment Record	001 KNEE PAIN
Number of Treatment Points	7

(b)(4)Proprietary Information



(b)(4)Proprietary Information





Copyright © 2012, *Theralase™ Inc.* 

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017







# TLC-1000 OPERATIONS MANUAL

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188

thera SE

TLC-1000 THERAPEUTIC MEDICAL LASER

(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11

thera SE

TLC-1000 THERAPEUTIC MEDICAL LASER

(b)(4)Proprietary Information	

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11

Page 2



## **1.** Safety Information – Please Review Before Using the Theralase System

The Theralase<sup>™</sup> TLC-1000 therapeutic laser is a Class 3B medical laser and as such is a potential hazard for direct and indirect viewing of the laser light.

The Theralase<sup>™</sup> TLC-1000 system as mentioned earlier is composed of three main components; the TLC-1000 Laser Controller, the Power Supply and the TLC-900 Multiple Laser Probe. Safety and safe use of the product are indicated both visually by the use of labels as well as equipment features. These will be detailed below. **IT IS IMPORTANT TO BE WELL AWARE OF ALL THESE SAFETY DETAILS PRIOR TO USING THE EQUIPMENT.** 

(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11

theralASE

b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11





**Operational Safety** 

Caution: Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure.

b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



b)(4)Proprietary Information

**TLC-1000 THERAPEUTIC MEDICAL LASER** 

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



(b)(4)Proprietary Information



## 2. THERALASE TLC-1000 SERIES PRODUCT INFORMATION

b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



## The Laser Controller (TLC-1000)

b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 



b)(4)Proprietary Information

**TLC-1000 THERAPEUTIC MEDICAL LASER** 

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 



#### CAUTION:

- Use only the TLC-154 power supply (+24 VDC, 2.7 A) to power the TLC-1000 laser controller. This is the only power supply approved to be used with the TLC-1000 laser controller.
- Outside of North America, the proper international cord set must be specified and purchased from Theralase <sup>™</sup> Inc. or purchased separately for the country of preference.



## The Multiple Laser Probe (TLC-900)

(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

thera SE

(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11




b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11





Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



b)(4)Proprietary Information



Hard case for TLC-1000 system (part number: TLC-1200)

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11

theralASE

b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



(b)(4)Proprietary Information

**TLC-1000 THERAPEUTIC MEDICAL LASER** 

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11

theral

b)(4)Proprietary Information

CDPH on 1 12 2017

+ 2016 6074

Figure 20. Outline of now the ILC-1000 system comes together

1. Unpack the contents of the case and identify each of the pieces as outlined in the figure above.



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11

theraldse

(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



(b)(4)Proprietary Information

thera SE

(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 



# 4. Operation Procedure (Knee Treatment Protocol)

(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 



	02 Options	5
(b)(4)	Proprietary	Informatior



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



b)(4)Proprietary Information

TLC-1000 THERAPEUTIC MEDICAL LASER

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017



TLC-1000 THERAPEUTIC MEDICAL LASER



(b)(4)Proprietary Information

TLC-1000 THERAPEUTIC MEDICAL LASER

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11

Page 35



b)(4)Proprietary Information



(b)(4)Proprietary Information



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 



b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11

theral

b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11

theral SE



Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11

Page 42



Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



# 7.Electro Magnetic Compatibility Information

In accordance with clause # 5.2.x of EN 60601-1-2.

Interconnecting Cables

b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



(b)(4)Proprietary Information

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11



b)(4)Proprietary Information

TLC-1000 THERAPEUTIC MEDICAL LASER

Copyright © 2012, *Theralase™ Inc.* 

July 10, 2012 - Rev. 1.11

Page 46

(b)(4)Proprietary Information



# b)(4)Proprietary Information

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017





# TLC-2000 OPERATIONS MANUAL

# TABLE OF CONTENTS

1.	SAFETY INFORMATION – PLEASE REVIEW BEFORE USING THE THERALASE SYSTEM	4
	Labelling and an Explanation of Symbols	4
	Location of the Labels and Symbols	4
	TLC-2000 Laser Controller Labels and Symbols	5
	Key Equipment Safety Features	6
	Hardware	
	Operational Safety	7
	Precautions	8
	Safety Warning	8
	Disclaimer Information	

2.	THERALASE TLC-1000 SERIES PRODUCT INFORMATION	
	Product Overview and Technical Description	
	The Power Supply (TLC-154)	
	The Laser Power Pack (TLC-2000)	
	Laser Output Probes (TLC-2001)	14
	Connecting/Disconnecting the Probes to the TLC-2000 Power Pack	
	Accessories	
	Safety Eyewear	
	Patient Safety	

START-UP INFORMATION	
The Power Pack	
The Probe	
The Tablet	
Set-up Procedure	23
	The Power Pack The Probe The Tablet

4.	SYSTEM OPERATION OVERVIEW	
	Operating the TLC-2000	
	How to Perform a Treatment	
	Stop Procedure	
	Overview of the Menu	
	System Options Tab	
	Practitioner Tab	
	Administrator Tab	
	System Operation in Detail	
	Common Operations:	
	Turning On/Off the Tablet:	
	Logon to the tablet:	
	Launching the TLC-2000 software:	
	Enter ng Pat ent Informat on	
	Treat ng the Pat ent	
	Ver fy ng the Laser Output	
	Advanced Operat ons	

How to Check Probe Usage Hours	
How to Ad ust Probe Buzzer Vo ume	
Adm n strator Funct ons	
Add/Remove users	
Set User Perm ss ons	
V ew Treatment records	
Warn ng Messages	

5		
	Inspect on and C ean ng	
	Week y	
	Month y	
	Prevent ve Ma ntenance Schedu e	
	How to check the opt ca output of each aser	
	Transport	
	Calibration/Repair of the Laser Verification Unit	
	Serv ce of Laser D odes	
	Product Warranty	
	Env ronmenta Cons derat ons	
	Thera ase Customer Serv ce Contact Informat on	

# 1. Safety Information – Please Review Before Using the Theralase System

The Theralase<sup>™</sup> TLC-2000 therapeutic laser is a Class 3B medical laser and as such is a potential hazard for direct and indirect viewing of the laser light.

The Theralase TLC-2000 system as mentioned earlier is composed of four main components; the tablet computer, the Power Supply, the Power Pack and the Laser Probes. The Laser Probes are the applied parts of the system. Safety and safe use of the product are indicated both visually by the use of labels as well as equipment features. These will be detailed below.

### NOTE: IT IS IMPORTANT TO BE WELL AWARE OF ALL THESE SAFETY DETAILS PRIOR TO USING THE EQUIPMENT.

Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

# Labeling and an Explanation of Symbols







TLC-2000 Laser Probe



TLC-2003 Power Supp y

F gure Locat on of the safety symbo s on the Thera ase TLC-2000 System

(b)(4)Proprietary Information

### **KEY EQUIPMENT SAFETY FEATURES**

### HARDWARE

In addition to the safety signs placed on each of the elements of the Theralase<sup>™</sup> system, key safety features are built into the TLC-2000 Power Pack and Tablet software as well. These are illustrated below.



Figure 3. Safety Features

**1. Emergency Stop:** This button shuts down the entire laser system immediately when pressed.

**2. ON/OFF Key Switch:** This is the main on / off switch for the laser system. When not in use, the keys can be placed in a safe location to prevent unauthorized use of the equipment.

**3. Remote Interlock Connector**: The Remote Interlock Connector is a connector that may be removed for disconnecting the mains power.

**4. Laser Probe:** The Laser Probes are the applied parts of the system. The Probe cannot operate the laser diodes unless a laser protocol has been received from the tablet software.

**5. Tablet Software:** Even if the Key Switch is left in the system with the Remote Interlock Connector in place and the Laser Probe(s) turned on, there is yet another security feature in place. The Laser Probe(s) cannot emit laser light unless treatment data has been received from the Theralase software on the tablet. To login into the operating system the Theralase TLC-2000 application software requires a username and password to operate.

### (b)(4)Proprietary Information

### **OPERATIONAL SAFETY**

Caution: Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure.

(b)(4)Proprietary Information

# SAFETY WARNING

- The TLC-2002 Laser Power Pack is designed to operate with a TLC-2001 Laser Probe only.
- The Theralase<sup>™</sup> TLC-2000 laser and all attachments are not designed to be protected from water intrusion and should never be used wet or in the accompaniment of water. If the laser system becomes wet or immersed in water, thoroughly dry the laser system and allow 48 to 72 hours in a warm environment to fully dry prior to use. If problems persist, contact Theralase<sup>™</sup> Inc. for service.
- The Theralase<sup>™</sup> TLC-2000 laser and all its attachments should not be used in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide.
- The Theralase<sup>™</sup> TLC-2000 laser and all attachments are not suitable for use in the presence of strong electromagnetic field or other electrical interference.
- The Theralase<sup>™</sup> TLC-2000 laser and all attachments are designed as an integral system and use of parts other than the recognized accessories or parts can degrade the performance of the Theralase<sup>™</sup> TLC-2000 laser and may result in hazardous radiation exposure.

WARNING: No Modification of This Equipment is Allowed WARNING: To Avoid Risk of Electrical Shock, This Equipment Must Only Be Connected to a Supply Mains With Protective Earth.

### DISCLAIMER INFORMATION

The TLC-2000 therapeutic laser system is a medical device and is sold to healthcare practitioners as a tool for the rehabilitation of their patients. The TLC-2000 therapeutic laser system and attachments are utilized to assist in the alleviation of their patient's symptoms. The TLC-2000 therapeutic laser system is not promoted as a cure for any patient ailments.

The information contained in this manual is for reference only and does not preclude the practitioner from using proper judgment in the rehabilitation of their patients.

It is the practitioner's sole responsibility to decide if laser therapy is indicated for their specific patient's rehabilitation.

The treatment protocols detailed in this manual are for reference only and are not promoted as the only protocol available to the practitioner.

The Theralase TLC-2000 Therapeutic Medical Laser System is "Indicated for Adjunctive Use in the Temporary Relief of Pain Associated with Knee Disorders with Standard Chiropractic Practice." (b)(4)Proprietary Information

(b)(4)Proprietary Information

# THE LASER POWER PACK (TLC-2002)

The Theralase<sup>™</sup> TLC-2000 laser system Power Pack (TLC-2002) is shown below. The Power Pack consists of an aluminum case with a plastic top and bottom.



Laser Probe Port



Unit ON Indicator



Laser Probe Power ON Indicator



Remote Interlock Connector



Emergency (E-STOP) Switch



Key Switch

Figure 5. The TLC-2000 system Power Pack (TLC-2002)

The Theralase<sup>™</sup> TLC-2000 laser system and all attachments are rated as Class 1 ME equipment, type BF and are equipped with Class 3B lasers.

Details regarding the safe use of the TLC-2000 system components are detailed in section 2 and operation of the unit in sections 3 and 4.

The top of the unit has an emergency STOP button (red button with white arrows), a key switch for switching the unit off and being able to lock the unit from unauthorized use, an interlock switch and a green ON LED. A further detail of these items as well as their operation is detailed in section 2.

The following are the detachable parts for the TLC-2000:

Power Supply:TLC-2003 Power Supply (+24 VDC, 4.2 A)Laser Probe:TLC-2001 Laser Probe

## WARNING: PLEASE DISCONNECT POWER SUPPLY BEFORE SERVICING

The Power Supply (TLC-2003)

The TLC-2003 is a universal power supply – the input voltage range is 100VAC to 240VAC, 47 to 63 Hz, meaning that it can be safely used anywhere in the world. The TLC-2003 power supply is able to power the TLC-2002 Power Pack.

Plugs in to main AC line via power cord (TLC-155)

24VDC 4A

Power output: Voltage: +24 VDC Current: 4.0 A



Power supply TLC-2003 plugs in to the Power Pack - TLC-2002

Figure 7. TLC-2003 Power Supply and associated Line Cord

With the proper cord set the TLC-2003 can be plugged in anywhere in the world. The cord set (TLC-155) supplied with the Theralase™ TLC-2003 power supply is a hospital grade 115 to 120 VAC, 50 to 60 Hz approved cord set. This is built to a North American standard and was chosen to work in Canada and the United States.

Some of the key features of the TLC-2003 Power Supply include

100-240 VAC Universal Input In Case IEC320 Single Output up to 101 Watts Outputs Regulated with Low Ripple Meets Safety Agency Requirements Complies with EMC/EMI Regulations CE Compliant Impact Resistant Polycarbonate Enclosure

### CAUTION:

• Use only the TLC-2003 Power Supply (+24 VDC, 4.2 A) to power the TLC-2002 Power Pack. This is the only Power Supply approved to be used with the TLC-2000 laser system.

• Outside of North America, the proper international cord set must be specified and purchased from Theralase<sup>®</sup> Inc. or purchased separately for the country of preference.

b)(4)Proprietary Information

### CAUTION

The Laser Probe Cover should be placed over the lens area when not in use.

### Connecting/Disconnecting the Laser Probe (TLC-2001) to the Power Pack (TLC-2002)

The Laser Probe features a very simple medical grade push-pull cable connector. This consists of a cable with a plug assembly on each end and a printed circuit board-mount receptacle assembly in the Power Pack and the Laser Probe.

Assurance of proper mating is provided by keying. The key(s) of the cable plug assembly are designated by an arrow on the plug and must align with the key of the receptacle assembly.




Figure 11. Probe Assembly illustrating the key

The cable plug assembly must be properly aligned with the receptacle assembly, and then inserted until the latches engage the receptacle assembly. There should be no need to force the connector into the receptacle. To ensure proper mating, there must be an audible "click". Refer to Figure 12.



Figure 12. Mating the Connectors - Cable Plug Assembly Key is aligned with Receptacle key

The connectors must be unmated by first gently sliding the housing subassembly away from the receptacle assembly – see Figure 13 and then pulling the entire connector out. PLEASE hold the main body and DO NOT pull the cable.




Figure 13. Un-mating the connected assembly

# ACCESSORIES

Accompanying the TLC-2000 system, in addition to the above main items, will be two sets of laser safety goggles (TLC-901).

### Safety Eyewear:


Two sets of laser safety goggles (part #: TLC-901) are supplied by Theralase as standard accessories. These laser safety goggles are provided for viewing of diffused laser light only and are not indicated for direct viewing of the laser output.

It is highly recommended that both the practitioner and the patient wear laser safety eyewear during laser therapy treatment.



Figure 14. Theralase supplied Laser Safety Eyewear

Optical Density ("O.D.") is the scientific term to describe a material's ability to block light of a certain wavelength. It is a logarithmic formula, where the higher the O.D. value, the better the material blocks light. For example, an O.D. of 1 means the material will allow 10% of the light penetrate the eyewear at a particular wavelength.



b)(4)Proprietary Information

Therefore, using the Theralase<sup>™</sup> approved laser glasses may result in the following levels of light being transmitted to the eye (assuming of course that the laser diode is directly in contact with the laser glasses and pointed at the eyes on the individual's face):

b)(4)Proprietary Information

### **CAUTION:**

Use only Theralase<sup>™</sup> Inc. supplied laser safety eyewear or laser safety eyewear that meets or surpasses the above noted optical densities at the above-specified wavelengths

NEVER POINT A LASER TOWARDS A PERSON'S EYES UNDER ANY CIRCUMSTANCES



Figure 14. Hard case for TLC-2000 system (part number: TLC-1200)

(b)(4)Proprietary Information	

# 3. START-UP INFORMATION

A detailed description of the individual system components

### The Power Pack (TLC-2002)

The TLC-2002 Power Pack provides power to the TLC-2001 Laser Probe.



Figure 16. A detailed look at the Power Pack

**1. Remote Interlock Connector** The Remote Interlock Connector is capable of disconnecting the TLC-2001 Laser Probe from main power when removed. To remove, press the button on the connector and pull the Remote Interlock Connector straight out. Main power is now disabled and no action will restore it until the Remote Interlock Connector is replaced. To replace, insert straight into the matching receptacle.

2. Emergency Stop This button shuts down the entire laser system immediately when pressed, disabling mains power.

# NOTE: Depressing the Red Emergency Stop Button immediately disables the laser system. Operation may be resumed by rotating the Red Emergency Stop Button 1/4 of a turn to the right and allowing it to "pop out" again.

**3. Key Switch** This is the main on / off switch for the laser system. When the key is turned to the left, the system is off, disabling mains power. When the key is turned ¼ turn to the right, the system is on. The key can only be removed when in the OFF position.



TLC-2003 Power Supply input



Connection to the TLC-2001

**4. Power Supply Input** Input for black DC adapter cord from the TLC-2003 Power Supply. When the DC adapter is connected the TLC-2000 Power Pack may be operated directly from an AC wall outlet.

Caution: Use only the TLC-2003 Power Supply for the Theralase<sup>™</sup> system.

5. Probe Power Port Power supply port for the TLC-2001 Laser Probe (High Density 6 pin medical grade connector).

**6. Power IN LED** This LED is illuminated only when the external power supply is connected and the key switch is in the ON position, the Remote Interlock Connector is installed and the Emergency Stop is not activated.

**7. Power OUT LED** This LED is illuminated only when both of the internal power supply voltages to the probes are good.

Note: The TLC-2001 Laser Probe will not function if this LED is OFF.

### The PROBE

The Theralase TLC-2001 Laser Probe is shown below. A feature of the Laser Probe is the use of a heatsink to keep the probe cool during use. There are three different strap sizes provided with the Laser Probes to facilitate holding the probe on different aspects of the knee.



Figure 16. A detailed look at the Laser Probe (TLC-2001)

**1. Display** A 2.4" Liquid Crystal Display **("LCD")** is used to display the status of the Laser Probe (idle/armed, laser operation( on/off ), BlueTooth<sup>®</sup> quality, laser treatment time, laser power, and Laser Probe error messages.

**2. Display Window** The display window is made from optically transparent acrylic plastic and can be cleaned with isopropyl alcohol.

**3. Power On/Off Switch** Slide this switch forward to turn on the Laser Probe. Slide the switch back to turn off the Laser Probe.



Caution: NEVER look directly into this area with the lasers on.

**6. Laser Output Window** The laser output window is an optically transparent acrylic plastic mounted to a ring. The window is removable by turning the ring left until it stops and lifting the ring straight off the probe. The window is coated to be scratch resistant. It keeps the lasers clean and can be wiped with isopropyl alcohol. Windows that become scratched should be replaced. Replacement windows can be ordered directly from Theralase.

# THE TABLET (TLC-2004)

The Tablet comes with the Theralase<sup>®</sup> software already installed and the keyboard and touch screen together. The touch screen can be separated from the keyboard if desired. There is a pop up touch screen keyboard; however, it may be easier to input data using the attached keyboard. With the keyboard attached, selections can be made using the cursor or the touch screen.

When delivered, the Tablet has an administrator username and password. It is recommended that the system be configured for new passwords and user names. Refer to the section Setting Username and Passwords in the Administrator Menu section for this procedure.

The Tablet is shipped with the manufacturer's user manual and software license card. Refer to those documents for additional information.



Figure 16. A detailed look at the tablet

- **1. Power On** Hold this button down for two (2) seconds until the touch screen turns on.
- **2. Charging Port** Plug the micro USB connector from the charger into this port to power the tablet from AC voltage.
- 3. USB Port Plug the Blue Tooth<sup>®</sup> USB stick into this port.
- **4. Release Button** Press this button to separate the keyboard and the touch screen.

**5. Touch Pad** Use this pad to control the cursor and make selections by moving your finger over the surface of the touch pad. Press the lower left corner of the touch pad to make a selection.

(b)(4)Proprietary Information

b)(4)Proprietary Information

### **REMINDER:**

The Theralase software must be operating with patient data in the database **BEFORE** the Laser Probe will emit laser light.

# **4. SYSTEM OPERATION OVERVIEW**

b)(4)Proprietary Information

### **Stop Procedure**

The Laser Probe may be turned off during operations by either of the following means:

- 1. Slide the switch on the Laser Probe to the OFF position.
- 2. Press the STOP button on the Tablet software
- 3. Press the Red Emergency Stop Button to stop Laser Probe operation and disable the mains power to the TLC-2000 laser system. Laser operation can only be re-started by turning the E-STOP switch a quarter turn clockwise to allow it to "pop" out.
- 4. Turn the on / off power key switch to the left ¼ turn to stop laser operation and to disable the mains power to the TLC-2000 laser system.

The Key Switch or Remote Interlock Connector may now be removed if desired to disable the laser system from unauthorized operation.

**Note:** Even if the Key Switch is left in the system and the Remote Interlock Connector is in place, the security password must be known to enable laser operation; thereby, providing a 3 step security protocol to prohibit laser use by unauthorized personnel.

### **Overview of the Menu**



b)(4)Proprietary

SYSTEM OPTIONS TAB

(b)(4)Proprietary Information

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188

### Turning the Tablet OFF

To turn the tablet OFF, swipe the screen in a Right to Left motion from the right hand edge of the screen. A vertical menu will be displayed. (To close the vertical menu, swipe the screen in a Left to Right motion.) Touch the Settings Icon to display the settings menu. Touch the power icon and select SHUT DOWN from the pop up window.

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188

QUESTIONS? CONTACT FDAVCDRH/OCE/DID at CDRH-FOISTATUS.GOV@IUA.INIS.gov of 301-790-0100

	Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-12-2017	
(b)(4)Proprietary Information		

(b)(4)Proprietary Information		

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS.GOV@fda.hhs.gov or 301-796-8188

### WARNING MESSAGES

b)(4)Proprietary Information
b)(4)Proprietary Information

### 5. Maintenance and Service

This section deals with maintenance and serviceability of the Theralase<sup>™</sup> system by the end user. If further detailed service of your laser system is required, contact Theralase Inc directly for assistance. Detailed contact information is supplied at the end of this section.

### **Inspection and Cleaning**

#### PRE AND POST TREATMENT

The probe should be cleaned with a 70% solution of isopropyl alcohol before and after treatments.

### WEEKLY

Follow the stop procedure in the basic operation section of this manual. Remove the key from the power pack. Disconnect the Laser Probe and power input to the Power Pack. Carefully clean the TLC-2001 Laser Probe with isopropyl alcohol swabs being careful not to immerse the Laser Probe head. Remove the lense cap and immerse a Q-tip into isopropyl alcohol and rotate in the laser aperture carefully to avoid damaging the laser diodes.

### MONTHLY

Follow the stop procedure in the basic operation section of this manual. Remove the key from the power pack. Disconnect the Laser Probe and power input to the power pack. Inspect the laser diodes for cracks or chips in the glass. A visible chip in the glass, a crack in the glass or a missing glass mandates replacement of the laser diode. Contact Theralase immediately for replacement of the defective / damaged laser diode.



### Transport

Remove the key from the main power switch and store in a safe location. Disassemble the laser system, carefully packing all components and taking special care with the laser probe assembly. In packing and unpacking the apparatus, be careful not to twist or kink the probe cable too sharply to prevent damage to the probe cable assembly. Ensure proper storage of the Laser Probe to avoid laser diode damage. If you damage your laser diodes, please notify Theralase<sup>™</sup> Inc. imme-diately for service / replacement. (The laser system should be kept within 0 to 40°C, 500 hPa to 1068 hPa and 0 to 80% relative humidity during transport).

#### Caution:

Calibration must only be performed by Theralase<sup>™</sup> Inc. or by trained personnel possessing the proper knowledge and test equipment. Contact Theralase<sup>™</sup> Inc., if after testing, software indicates maintenance of the laser diodes is required.

### **Environmental Considerations**

This product contains electronic components and as such is considered e-waste and must be disposed of responsibly. Do not dispose of this product with waste destined for landfill or incineration. This product must be recycled according to local laws & regulations. At the end of this product service life, estimated to be 10 years, the product can be returned to Theralase for disposal.

Records Processed under FOIA Request 2016-6074. Released by CDRH (b)(4) Proprietary Information -

Client	Theralase Inc.
Product	TLC-2000
Standard(s)	(b)(4)Proprietary Information -

# Report Scope

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

### Summary

The results contained in this report relate only to the item(s) tested. (4)Proprietary Information - Testing

For testing dates, see "Testing Environmental Conditions and Dates".

Product     TLC-2000       Standard(s)     (b)(4)Proprietary Information -	Client	Theralase Inc.	Testing
Standard(s) (b)(4)Proprietary Information -	Product	TLC-2000	
	Standard(s)	(b)(4)Proprietary Information -	

### Test Results Summary

(b)(4)Proprietary Information - Testing	

Records Processed under FOIA Request 2016-6074. Released by CDRH on (b) (4) Proprietary Information -

Client	Theralase Inc.
Product	TLC-2000
Standard(s)	(b)(4)Proprietary Information -

 Records Processed under FOIA Request 2016-6074. Released by CDRI
 (b)(4)Proprietary Information - Testing

 Client
 Theralase Inc.

 Product
 TLC-2000

 Standard(s)
 (b)(4)Proprietary Information - Testing

# Applicable Standards, Specifications and Methods

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

(b)(4)Proprietary Information -	Testing	

Client	Theralase Inc.	Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information - Testing	

# **Document Revision Status**

Revision 0: April 20, 2015 First Release

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

# Definitions and Acronyms

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

# **Testing Facility**

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

### **Testing Environmental Conditions and Dates**

Following environmental conditions were recorded in the facility during time of testing

(b)(4)Proprietary Information	ation - Testing		

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

### **Test Results Section**

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

#### b)(4)Proprietary Information - Testing

Records Processed under FOIA Request 2016-6074. Released by CDRH c(b)(4)Proprietary

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on the second se

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH c(b)(4) Proprietary Information -

Client	Theralase Inc.	
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

b)(4)Proprietary Information - Testing

1

Client	Theralase Inc.	(b)(4)Proprietary Information Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on 101(4) Proprietary

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH o(b)(4) Proprietary Information -

Client	Theralase Inc.	Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on by 2 Proprietary

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

|--|

	Client	Theralase Inc.	Information - Testing
	Product	TLC-2000	
	Standard(s)	(b)(4)Proprietary Information -	
priet	ary Information -	Testing	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1 10 20 Proprietary

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	-
Standard(s)	(b)(4)Proprietary Information -	-

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on the Proprietary

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

# Test Equipment List

(b)(4)Proprietary Information - Testing	

Client	Theralase Inc.	(b)(4)Proprietary Information Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1-160447 Proprietary

Client	Theralase Inc.	Information - Tes
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	- Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

(b)(4)Proprietary Information -	Testing

Records Processed under FOIA Request 2016-6074. Released by CDRH on the relation -

Client	Theralase Inc.	Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

[	Client	Theralase Inc.	(b)(4)Proprietary Information Testing
[	Product	TLC-2000	-
	Standard(s)	(b)(4)Proprietary Information -	а -

Records Processed under FOIA Request 2016-6074. Released by CDRH or (b)(4)Proprietary Information

Client	Theralase Inc Testing	
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH o(b)(4) Proprietary Information -

Client	Theralase Inc.	Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	IEC/EN 60601-1-2:2007 - (Ed3.0)	

Records Processed under FOIA Request 2016-6074. Released by CDRH on 101(4)Proprietary

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on (b)(4) Proprietary

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on(c)(4)Proprietary

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1 (5) (4) Proprietary Information

Client	Theralase Inc.	- Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH (b)(4)?:oprietary Information -

Client	Theralase Inc.	Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH or 1/4 Proprietary

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	а 

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	а 

Records Processed under FOIA Request 2016-6074. Released by CDRH (b)(4) Proprietary Information -

Client	Theralase Inc.	lesting
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Informatio
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)		

Client	Theralase Inc.	Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on the Proprietary Information -

Client	Theralase Inc.	Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	-

Records Processed under FOIA Request 2016-6074. Released by CDRH on (1) (4) Proprietary Information

Client	Theralase Inc.	- Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH or (b) (4) Proprietary

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on (b)(4) Proprietary Information

Client	Theralase Inc.
Product	TLC-2000
Standard(s)	(b)(4)Proprietary Information -

Records Processed under FOIA Request 2016-6074. Released by CDRH on (b)(4) Proprietary Information -

Client	Theralase Inc.	Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Re	cords Processed under FOIA Request 2016-6074. Released by CDRH on	(b)(4)Proprietary
Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	- Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

	Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
	Product	TLC-2000	
	Standard(s)	(b)(4)Proprietary Information -	
b)(4)Pro	oprietary Informat	ion - Testing	

Records Processed under FOIA Request 2016-6074. Released by CDRH or (b)(4)Proprietary Information -

Client	Theralase Inc.	Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	a

Client	Theralase Inc.	(b)(4)Proprietary Informatio - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

10

Client	Theralase Inc.       (b)(4)Proprietary         Information - Testing
Product	TLC-2000
Standard(s)	(b)(4)Proprietary Information -

Records Processed under FOIA Request 2016-6074. Released by CDRH on the Vertice and Informatic

Client	Theralase Inc.	Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Product     TLC-2000       Standard(s)     (b)(4)Proprietary Information -	etary - Testing
Standard(s) (b)(4)Proprietary Information -	
Taction	

Client	Theralase Inc.	(b)(4)Proprietary Information Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on the 14 Proprietary

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

(b)(4)Proprietary Information - Testing	

Client	Theralase Inc.	Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Infor

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH or (b)(4). reprietary Information -

Client	Theralase Inc.	
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	
(4)Proprietary Inf	ormation - Testing	

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on (b)(4)Proprietary

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on (b)(4)Proprietary Informatio

Client	Theralase Inc.	- Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on (b)(4)Proprietary

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1 (b) (4) Proprietary

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information · Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

- 1

	Client	Theralase Inc Testing	etary mormation
	Product	TLC-2000	
	Standard(s)	(b)(4)Proprietary Information -	
(4)	Proprietary Inform	mation - Testing	

Client	Theralase Inc.	(b)(4)Proprietary Information Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Informa - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH or (b)(4) Proprietary Information -

Client	Theralase Inc.
Product	TLC-2000
Standard(s)	(b)(4)Proprietary Information -

Client	Theralase Inc.	(b)(4)Proprietary Information Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

© Global EMC Inc This test report shall not be reproduced except in full, without written approval of Global EMC Inc

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

- 1

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on the 14 30 Toprietary Information

Client	Theralase Inc.	- Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

(4)Proprietary Information - Testing		

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	

Records Processed under FOIA Request 2016-6074. Released by CDRH on 1 (b) (4) Proprietary

Client	Theralase Inc.	Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	en e

Client	Theralase Inc.	(b)(4)Proprietary Information - Testing
Product	TLC-2000	
Standard(s)	(b)(4)Proprietary Information -	