



# U.S. Department of Health & Human Services

---

**Food and Drug Administration**

## SAVE REQUEST

**USER:** (kml)  
**FOLDER:** K941507 - 303 pages  
**COMPANY:** SANESE MEDICAL CORP. (SANEMEDI)  
**PRODUCT:** LAPAROSCOPE, GYNECOLOGIC (AND ACCESSORIES) (HET)  
**SUMMARY:** Product: THEMO-FLO SYSTEM 3  
LAPAROSCOPIC/HYSTEROSCOPIC SYSTEM  
**DATE REQUESTED:** Oct 5, 2015  
**DATE PRINTED:** Oct 5, 2015  
**Note:** Printed



K941507

SECTION 510(K) SUMMARY OF SAFETY AND  
EFFECTIVENESS INFORMATION FOR SMC'S  
THERMO-FLO SYSTEM 3 LAPAROSCOPIC/  
HYSTEROSCOPIC IRRIGATION SYSTEM AND  
SANESE DISPOSABLE IRRIGATION/ASPIRATION  
PROBE AND TUBING

Pursuant to Section 12(a) of the Safe Medical Devices Act of 1990, 21 U.S.C. § 360c(i)(3)A and 21 C.F.R. § 807.92, the following summary of information concerning the safety and effectiveness of SMC's Thermo-Flo System 3 Laparoscopic/Hysteroscopic Irrigation System and Sanese Disposable Irrigation/Aspiration Probe and Tubing is hereby submitted.

1. Name and address of device manufacturer submitting 510(k) notification:  

Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, OH 43212  
(614) 291-8848
2. Regulatory correspondent.  

Christopher N. Sanese  
Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, OH 43212  
(614) 291-8848
3. Date summary was prepared: March 17, 1994
4. Names of Devices:
  - (i) Proprietary Names:
    - (a) Thermo-Flo System 3 Laparoscopic/Hysteroscopic Irrigation System\*
    - (b) Sanese Disposable Irrigation/Aspiration Probe and Tubing\*
  - (ii) Common, usual names:
    - (a) Laparoscopic irrigation pump system
    - (b) Disposable aspiration/irrigation probe and tubing.
  - (iii) Classification Names: Hysteroscope and Accessories and Gynecologic Laparoscope and Accessories, 21 C.F.R. §§ 880.1690 and 884.1720.
5. Predicate devices to which SMC is claiming substantial equivalence:

SMC claims substantial equivalence for purposes of section 510(k) only to:

A. Irrigation Pump System:

The Niagra™ High-Flow Irrigator System, 510(k) Number K924530 by Cabot Medical, 2021 Cabot Boulevard, West Langhorn, PA 19047 (800-523-6078).

B. Disposable Probe and Tubing:

(1) Probe:

Corson™ Disposable Suction Irrigation Probe by Cabot Medical, 2021 Cabot Boulevard, West Langhorn, PA 19047 (800-523-6078).

(2) Tubing:

Core Dynamics Inc. Disposable Irrigation/Suction Tubing Set (510(k) Number K910246) by Core Dynamics, Inc., 9951 Atlantic Boulevard, Jacksonville, Florida 3225 (904-727-0910).

6. Description of the devices and primary principles of operation:

A. Description of irrigation pump system.

This device consists of a fabric enclosure (pressure envelope) in which a bag of irrigation fluid is placed. The fabric enclosure contains an inflatable bladder that is pneumatically inflated to pressurize the fluid in the irrigation fluid bag. Inflation of the bladder is provided by an air compressor which has a maximum voltage leakage current of 100 microamps. This complies with Underwriters Laboratory Standard 544.

The pressure envelope also contains a heating element to assist in maintaining the temperature of preheated irrigation solution while being used in the pressure envelope. There is also a thermal regulator to control the temperature of the heating element. The heating element and the thermal regulator are powered by an Underwriter Laboratory approved 24 volt electrical power supply.

The air compressor and heating element are controlled by a microcontroller. These components are controlled by the micro-controller through activating solid state switches which control electrical power to the pump, pump solenoids and the heater.

The pressurized preheated fluid is then introduced into the endoscopic surgery site through a medical grade plastic tube and probe that are discussed below.

B. Description of disposable probe and tubing:

The probe consists of a medical grade plastic, activation valve and metal shaft that will deliver and aspirate fluids to and from an operative site.

The tubing consists of two 15 foot length medical grade surgical tubing that are approximately 10 millimeters in diameter which would be connected to the probe and the irrigation solution bag and the pump when used for irrigation or to the probe and the suction unit when used for aspiration. The probe and tubing are made of medical grade plastic and metal that have been determined by FDA to be safe and effective for this intended usage.

C. Primary Principles of Operation:

(i) Irrigation Pump System:

The primary principles of operation are pressurization through an electrically powered pneumatic compressor for pressurized irrigation fluid flow and thermal transfer through an electrically-powered thermo-conducted heat transfer unit to maintain preheated fluid.

(ii) Probe and Tubing:

With respect to the aspiration/irrigation use of the probe and tubing, the principle operation is vacuum and pressurization.

7. Statement of the purpose for which the device will be recommended:

These devices are intended to be used for irrigation and aspiration during adult and pediatric hysteroscopic or laparoscopic (endoscopic) surgery.

8. Biocompatibility and effectiveness of the device:


(i) Biocompatibility:

The only aspect of the device that comes in contact with the patient is the distal end of the probe. The probe and the probe tubing have been subject to biocompatibility tests and are safe and effective for such usage. The irrigation pump has a maximum voltage



leakage current of 100 microamps in accordance with U.L. Standard 544. Performance tests of the pump demonstrate that it is safe and effective for its intended use. Also, all components of the pump meet U.L. Standards.

Respectfully submitted,



Christopher N. Sanese  
Regulatory Correspondent  
Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, OH 43212  
(614) 291-8848

**SANESE MEDICAL**<sup>®</sup>  
CORPORATION

MANUFACTURER OF MEDICAL PRODUCTS

K941507/A<sup>3</sup>  
**RECEIVED** ADD TO FILE  
7 MAR 95 11 42  
FDA/CDRH/ODE/DMC

March 3, 1995

Ms. Donna-Bea Tillman  
FDA  
9200 Corporate Boulevard  
Rockville, Maryland 20850

RE:K941507 / A  
Thermo-flo System 3  
Dated: July 5, 1994  
Received: July 6, 1994


VIA FAX & FEDEX  
301-594-2339

Dear Donna,

Per our telephone conversation of today, I am forwarding the revised draft labeling you requested.

Thank you for your time and consideration.


Sincerely,




Christopher N. Sanese  
Regulatory Corespondent  
SANESE MEDICAL CORPORATION

ENC: 1  
cs/CNS

cc: John Jevicky, Esq.  
DINSMORE & SHOHL



# OPERATION MANUAL



**DRAFT COPY**

**INDICATION FOR USE:** INDICATED FOR IRRIGATION DURING LAPAROSCOPIC PROCEDURES.

**CONTRA INDICATIONS:** WHERE LAPAROSCOPIC SURGERY IS CONTRA INDICATED.

**FEDERAL (USA) LAW** restricts this device to sale by, or on the order of a physician. Recommended for Laparoscopic Irrigation in adult and pediatric surgery.

**THIS PRODUCT TO BE USED FOR LAPAROSCOPIC IRRIGATION. DO NOT USE FOR HYSTEROSCOPIC DISTENTION.**

**SANESE MEDICAL**  
CORPORATION



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

MAR - 6 1995

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mr. Christopher N. Sanese  
Regulatory Correspondent  
Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, Ohio 43212

Re: K941507/B  
Thermo-Flo System 3  
Laparoscopic/Hysteroscopic  
Irrigation System®  
Sanese Disposable  
Irrigation/Aspiration Probe  
and Tubing™  
Dated: January 13, 1995  
Received: January 17, 1995  
Regulatory class: II  
21 CFR §884.1720/Procode: 85 HET

Dear Mr. Sanese:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug and Cosmetic Act (Act). You may therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval) it may be subject to such additional controls. Existing major regulations affecting your device can be found in the code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, FDA will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, the FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you may have under the Electronic Product Radiation Control provisions, or other Federal Laws or Regulations.

This letter immediately will allow you to begin marketing your device as described in your 510(k) premarket notification. An FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and permits your device to proceed to the market, but it does not mean that FDA approves your device. Therefore, you may not promote or in any way represent your device or its labeling as being approved by FDA. If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), promotion, or advertising please contact the Office of Compliance, Promotion and Advertising Policy Staff (HFZ-302) at (301) 594-4639. Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll free number (800) 638-2041 or at (301) 443-6597.

Sincerely yours,

Lillian Yin, Ph.D.  
Director, Division of Reproductive,  
Abdominal, Ear, Nose, and Throat,  
and Radiological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

7a 1/24/95 RR

510(K) ROUTE SLIP

510(k) NUMBER K941507 PANEL OB DIVISION DRAER BRANCH OGDB

TRADE NAME THEMO-FLO SYSTEM 3 LAPAROSCOPIC/HYSTEROSCOPIC SYSTEM

COMMON NAME LAPAROSCOPIC IRRIGATION PUMP SYSTEM

PRODUCT CODE \_\_\_\_\_

APPLICANT SANESE MEDICAL CORPORATION

SHORT NAME SANEMEDI

CONTACT CHRISTOPHER N SANESE

DIVISION \_\_\_\_\_

ADDRESS 885 NORTHWEST BLVD.

COLUMBUS, OH 43212

PHONE NO. (614) 291-8848

FAX NO. ( ) - -

MANUFACTURER SANESE MEDICAL CORPORATION

REGISTRATION NO. \_\_\_\_\_

DATE ON SUBMISSION 16-MAR-94

DATE DUE TO 510(K) STAFF 12-JUN-94

DATE RECEIVED IN ODE 29-MAR-94

DATE DECISION DUE 27-JUN-94

DECISION \_\_\_\_\_

DECISION DATE MAR - 6 1995

SUPPLEMENTS	SUBMITTED	RECEIVED	DUE POS.	DUE	OUT
<u>S002</u>	<u>13-JAN-95</u>	<u>17-JAN-95</u>	<u>02-APR-95</u>	<u>17-APR-95</u>	
<u>S001</u>	<u>05-JUL-94</u>	<u>06-JUL-94</u>	<u>19-SEP-94</u>	<u>04-OCT-94</u>	<u>17-NOV-94</u>

CORRESPONDENCE	SENT	DUE BACK	
<u>C001</u>	<u>06-JUN-94</u>	<u>06-JUL-94</u>	<u>HOLD LETTER</u>
<u>C002</u>	<u>17-NOV-94</u>	<u>17-FEB-95</u>	<u>HOLD LETTER</u>

OTHER SUBMISSIONS	SUBMITTED	RECEIVED	DUE POS	DUE	CUT
<u>ADD-TO-FILE</u>	<u>10-AUG-94</u>	<u>15-AUG-94</u>			
<u>ADD-TO-FILE</u>	<u>15-OCT-94</u>	<u>19-OCT-94</u>			

2

Handwritten signature and scribbles covering the bottom half of the page.



# Memorandum

Date March 3, 1995

From REVIEWER(S) - NAME(S) Donna-Bea Tillman

Subject 510(k) NOTIFICATION K941507/S<sup>2</sup>

To THE RECORD

It is my recommendation that the subject 510(k) Notification:

- (A) Is substantially equivalent to marketed devices.
- (B) Requires premarket approval. NOT substantially equivalent to marketed devices.
- (C) Requires more data.
- (D) Other (e.g., exempt by regulation, not a device, duplicate, etc.)

Additional Comments:

Is this device subject to Postmarket Surveillance? Yes  No

This 510(k) contains: (check appropriate box(es))

- A 510(k) summary of safety and effectiveness, or
- A 510(k) statement that safety and effectiveness information will be made available
- The required certification and summary for class III devices

The submitter requests under 21 CFR 807.95:\*

- No Confidentiality
- Confidentiality for 90 days
- Continued Confidentiality exceeding 90 days

Predicate Product Code w/panel and class:

85 HET Class II  
884.1720

Additional Product Code(s) w/Panel (optional):

REVIEW: Cheri M. Pollard  
(BRANCH CHIEF)


06DB | 3/6/95  
(BRANCH CODE) (DATE)

FINAL REVIEW: P. S. [Signature] / a Lin  
(DIVISION DIRECTOR)

3/6/95  
(DATE)

\*DOES NOT APPLY TO ANY "SE" DECISIONS

# 510(k) "SUBSTANTIAL EQUIVALENCE" DECISION-MAKING PROCESS (DETAILED)



- 510(k) submissions compare new devices to marketed devices. FDA requests additional information if the relationship between marketed and "predicate" (pre-Amendments or reclassified post-Amendments) devices is unclear.
- This decision is normally based on descriptive information alone, but limited testing information is sometimes required.
- Data may be in the 510(k), other 510(k), the Center's classification files, or the literature.

# OGDB

**K941507/B**

**Reviewer:** Donna-Bea Tillman  
Biomedical Engineer

**Division/Branch:** DRAERD/ADOU/OGDB  
(HFZ-470)

**Proprietary Trade Name:** Thermo-Flow System 3  
**Common Name:** Laparoscopic Irrigation and Aspiration

**Product to which compared:** Cabot Niagara High-Flow Irrigator (K924530)

**Applicant:** Sanese Medical Corporation  
885 Northwest Blvd  
Columbus, OH 43212

**Contact:** Christopher N. Sanese, Regulatory Correspondent  
**Phone:** (614) 291-8848

## DEVICE DESCRIPTION

1. *Intended Use:*

The Thermo-Flow System 3 is intended for irrigation/aspiration during laparoscopic surgery (Note: In a October 15, 1994, amendment, the firm withdrew its claims for hysteroscopy).

2. *Physical Description:*

	YES	NO
• Is the device life-supporting or life sustaining?	—	<u>X</u>
• Is the device implanted (short-term or long-term)?	—	<u>X</u>
• Does the device design use software?	<u>X</u>	—
• Is the device sterile?	—	<u>X</u>
• Is the device single use?	—	<u>X</u>
• Is the device home use?	—	<u>X</u>
• Is the device for prescription?	<u>X</u>	—
• Does the device contain a drug or biological product as a component?	—	<u>X</u>
• Is this device a kit?	—	<u>X</u>



(b)(4) Trade Secret Process - Product Specifications



## REVIEW ANALYSIS

(b)(4) Trade Secret Process - Review Analysis



4

(b)(4) Trade Secret Process - Review Analysis




(b)(4) Trade Secret Process - Review Analysis



8

(b)(4) Trade Secret Process - Review Analysis



9

## Substantial Equivalence (SE) Decision Making Documentation

	YES	NO	
1. IS PRODUCT A DEVICE?	<u>X</u>	—	IF NO, STOP
2. DEVICE SUBJECT TO 510(k)?	<u>X</u>	—	IF NO, STOP
3. SAME INDICATION STATEMENT?	<u>X</u>	—	IF YES, GOTO 5
4. DO DIFFERENCES ALTER THE EFFECT OR RAISE NEW ISSUES OF SAFETY OR EFFECTIVENESS?	—	—	IF YES, STOP->NSE
5. SAME TECHNOLOGICAL CHARACTERISTICS?	—	<u>X</u>	IF YES, GOTO 7
6. COULD THE NEW CHARACTERISTICS AFFECT SAFETY OR EFFECTIVENESS?	—	—	IF YES, GOTO 8
7. DESCRIPTIVE CHARACTERISTICS PRECISE ENOUGH?	—	<u>X</u>	IF YES, STOP->SE
8. NEW <u>TYPES</u> OF SAFETY AND EFFECTIVENESS QUESTIONS?	—	<u>X</u>	IF YES, STOP->NSE
9. ACCEPTED SCIENTIFIC METHODS EXIST?	<u>X</u>	—	IF NO, STOP->NSE
10. PERFORMANCE DATA AVAILABLE?	<u>X</u>	—	IF NO, REQ. DATA
11. DATA DEMONSTRATE EQUIVALENCE	—	<u>X</u>	IF YES, STOP->SE

10

# REVIEWER RECOMMENDATION

## Substantially Equivalent

ProCode: 85 HET  
Class: II  
CFR #: 21 CFR §884.1720

Donna Bea Tillman 3/3/95  
Donna-Bea Tillman, Ph.D Date

Concur

W. Wood  
3/6/95

Colin M. Pollard 3/6/95  
Colin M. Pollard Date  
Chief, Ob/Gyn Devices Branch

/✓ Concur  
/ / Do not concur.

Comments:

//

**SANESE MEDICAL®**  
CORPORATION

MANUFACTURER OF MEDICAL PRODUCTS

**DATE: March 3, 1995**  
**TO: Ms. Donna-Bea Tillman**  
**FIRM: FDA**  
**TELECOPIER NO: 301-594-2339**  
**FROM: Christopher N. Sanese**

NO. OF PAGES (INCLUDING COVER SHEET) TRANSMITTED: ~~X~~ 3

**IF THERE ARE ANY PROBLEMS IN RECEIVING THIS TRANSMISSION,  
PLEASE CALL (614)291-8848 IMMEDIATELY**

**THANK YOU,  
TELECOPIER OPERATOR**

SENT BY: \_\_\_\_\_  
TIME SENT: \_\_\_\_\_  
COMPLETED & LOGGED  
NOTES & COMMENTS:

*This message is intended only for the use of the individual or entity to which it is addressed and may contain information that is privileged, confidential and exempt from disclosure under applicable law. If the reader of this notice is not the intended recipient, you are hereby notified that any dissemination, distribution or copying of this communication is strictly prohibited. If you have received this communication in error, please notify us immediately by telephone and return these papers to us at the address shown above via first class mail.*

- \_\_\_\_\_ *No original documents will follow unless requested*
- \_\_\_\_\_ *Original document will follow by courier*
- \_\_\_\_\_ *Original documents will follow by mail*

12

**SANESE MEDICAL**<sup>®</sup>  
CORPORATION

MANUFACTURER OF MEDICAL PRODUCTS

March 3, 1995

Ms. Donna-Bea Tillman  
FDA  
9200 Corporate Boulevard  
Rockville, Maryland 20850

RE:K941507 / A  
Thermo-flo System 3  
Dated: July 5, 1994  
Received: July 6, 1994


VIA FAX & FEDEX  
301-594-2339

Dear Donna,

Per our telephone conversation of today, I am forwarding the revised draft labeling you requested.

Thank you for your time and consideration.

Sincerely,



Christopher N. Sancsc  
Regulatory Correspondent  
SANESE MEDICAL CORPORATION

ENC: 1  
cs/CNS


cc: John Jevicky, Esq.  
DINSMORE & SHOHL

13



# THERMO-FLO SYSTEM 3

## OPERATION MANUAL



**DRAFT COPY**

**INDICATION FOR USE:** INDICATED FOR IRRIGATION DURING LAPAROSCOPIC PROCEDURES.

**CONTRA INDICATIONS:** WHERE LAPAROSCOPIC SURGERY IS CONTRA INDICATED.

**FEDERAL (USA) LAW restricts this device to sale by, or on the order of a physician. Recommended for Laparoscopic Irrigation in adult and pediatric surgery.**

**THIS PRODUCT TO BE USED FOR LAPAROSCOPIC IRRIGATION. DO NOT USE FOR HYSTEROSCOPIC DISTENTION.**

74

**SANESE MEDICAL**

**SANESE MEDICAL**<sup>®</sup>  
CORPORATION

MANUFACTURER OF MEDICAL PRODUCTS

March 2, 1995

Ms. Donna-Bea Tillman  
FDA  
9200 Corporate Boulevard  
Rockville, Maryland 20850

RE:K941507 / A

Thermo-flo System 3

VIA FAX & FEDEX

Dated: July 5, 1994

Received: July 6


301-594-2339

Dear Donna,

Per our telephone conversation of February 24, 1995 I am forwarding the draft labeling copy you requested; A) Operators Manual B) Reusable probe label.

Thank you for your time and consideration.

Sincerely,



Christopher N. Sanese  
Regulatory Correspondent  
SANESE MEDICAL CORPORATION

ENC: 9

cs/CNS


cc: John Jevicky, Esq.

DINSMORE & SHOHL

---

**THERMO-FLO**  
SYSTEM 3

**OPERATION MANUAL**



**DRAFT COPY**

**SANESE MEDICAL**  
CORPORATION

16

**THERMO-FLO**  
SYSTEMS

**TABLE OF CONTENTS**

**1.0 Product Description**

**2.0 Safety**

**2.1 Symbol Explanations**

**2.2 Safety Precautions**

**2.3 Safety Features**

**3.0 Specifications**

**4.0 Operating The System**

**4.1 Description of Controls**

**4.2 Description of Operation**

**5.0 Usage/Care/Warranty**

**SANESE MEDICAL**  
CORPORATION

17

# 1.0 Product Description


The Thermo-Flo System 3 was developed to eliminate the major problems faced everyday by the surgeons and O.R. unit Managers. These are cross contamination, inadvertent patient hypothermia, inadequate irrigation solution capacity, ease of set up, tube sets and disposable probes that leak, noise and physical space. Sanese Medical answers these problems with the Thermo-Flo System 3.

1) **CROSS CONTAMINATION** occurs when non sterile fluid enters the patient's abdominal cavity through the irrigation line. The Thermo - Flo System 3 eliminates cross contamination.

2) **INADVERTENT PATIENT HYPOTHERMIA** is the result of the combination of gas and fluids ranging in temperature from 65 to 72 degrees fahrenheit being introduced into the patient. The Thermo - Flo System can maintain the solution temperature at 100 degrees fahrenheit. Though this system alone cannot prevent hypothermia, it will dramatically effect the reduction of the core temperature due to existing conditions of surgery.

### 3) SOLUTION CAPACITY / EASE OF SET UP

The Thermo - Flo System 3 meets the liter needs of the today's and tomorrow's laparoscopic procedures. This amount provides less disruption of the procedure, by not replacing the solution, cap and tube inserts that promotes the inadvertent contamination by the person replacing the solution as often.



- KEY ADVANTAGES OF THE THERMO-FLO SYSTEM 3**
- \* Delivers a strong stream of irrigation fluid
  - \* No capital equipment necessary
  - \* Choice of disposable or reusable probe
  - \* Simplified O.R. set up / LV.
  - \* Avoids Equipment purchase
  - \* Quiet operation
  - \* 24 hour, 7 days a week service available

## THERMO-FLOW SYSTEM 3 PRODUCT LINE

- 1) The Thermo-Flo System 3 Irrigation Pump.....TFS3 one per case
- 2) TFS3-3WP Tube Set with Disposable Probe.....TFS3-3WP one per case
- 3) TFS3-3WOP Tube Set with out Probe .....TFS3-3WOP one per case
- 4) TFS3-N-3WP Tube Set with Disposable Probe.....TFS3-N-3WP 10 per case
- 5) TFS3-n-3WP Tube set with out Probe.....TFS3-N-3WOP 10 per case
- 6) TFS3- 3RPR Reusable Probe Set .....TFS3-N-3RPR one per case

All tube sets are 15' in length .

**SANESE MEDICAL**  
CORPORATION  
390 NORTHWEST BLDG. COLUMBUS, OHIO 43212 (614) 291-4242 FAX (614) 291-1222

78

## 2.0 SAFETY

### 2.1 SYMBOL EXPLANATIONS



**READ THE OPERATIONS MANUAL BEFORE OPERATING THIS UNIT.  
FOR SALE BY, OR ON THE ORDER OF A HEALTH CARE PROFESSIONAL.**



**THIS UNIT CAN NOT BE OPERATED IN THE PRESENCE  
OF AN EXPLOSIVE ANAESTHETIC MIXTURE WITH AIR, OXYGEN  
OR NITROUS OXIDE.**



**Only qualified personnel should perform repair to this unit.  
This unit requires a power input of - 110 VAC 60 HZ 1.0 AMP  
This device is equipped with a "Hospital Grade" plug. Grounding  
reliability can only be achieved when the equipment is connected  
to an equivalent receptacle marked "Hospital Only" or "Hospital Grade".**

### 2.2 SAFETY PRECAUTIONS

Use only with TFS 3 irrigation tube sets. These sets are equipped with lock and pierce components for the safe operation of the device. Never use parts from another manufacturer's unit even though they appear to fit.

Never service device while it is plugged into an outlet.

Never service device near water.

This device is designed to be operated from a power indicated on the marking label.

Do not overload wall outlets or use extension cords.

Do not allow anything to rest on the power cord or locate this product where the cord will be abused by persons or equipment moving on or near the cord.

Unplug the unit from the wall and refer servicing to a qualified personnel under the following conditions:

- Damaged or frayed power cord
- If product has been dropped or case is damaged
- If a marked change in product performance is noted

### 2.3 SAFETY FEATURES

- OVER PRESSURE -** In the event over pressure occurs the unit cut power to the pump and air supply system.
- OVER TEMPERATURE -** In the event an over temperature condition occurs an independent hardware/software circuit will cut power to the heater.
- LEAKAGE CURRENT-** The leakage current of the pressure enclosure is less than 10 micro amps.
- WATCHDOG TIMER-** If the micro controller fails to execute code properly the watchdog timer will reset the Unit to the initial conditions.

19

### **3.0 SPECIFICATIONS Model TFS3 \***

**Power**

120 v ac 1 AMP or 220 V AC .5 AMP 50/60 Hz  
(95 to 250 V AC to 63 Hz)

**Weight**

7 lbs. (with out stand)

**Pressure Range**

100 mm Hg to 800 mm Hg + or - 10 mm Hg

**Heater Temperature**

The heater is preset to 100 F° + or - 2 F°  
Safety Shutdown Temperature- 105 F°.

**Maximum Leakage Current**

Control Unit - 100 micro amps. (per UL 544)

Pressure Enclosure - 10 micro amps. (per UL 544)

Maximum Ground Impedance - 0.1 ohm ( per UL 544)

Fuse : Slow Blow 2.0 amp. ( Located on the back of control unit)

**Environmental Operating Conditions**

40 F° to 150 F°, Humidity non-condensing

\* U. S. and International patent pending  
Specification subject to change without notice.

20

## 4.0 OPERATING THE SYSTEM

### 4.1 DESCRIPTION OF THE CONTROLS

**Power Switch:** The power switch is located on the front panel. It controls the power to unit.

**Pressure Adjust Switch:** Located to the right of the displays a three position rocker switch which is used to adjust the pressure. The center position is off. The top position (up arrow) and bottom position (down arrow) are momentary. Pressing the UP arrow increases the pressure. Pressing the bottom (down arrow) will decrease the pressure.

**Load/Run Switch:** This switch is a two position rocker. When the switch is set to the Load position the device will drop the air pressure in the Pressure Enclosure so the bag can be reloaded. When the switch is in the Run position the Pressure Enclosure will inflate to the set point on the display.

**Display :** The display is a 16 character single line LCD.

**Low Fluid LED:** There is a panel light on the front panel that indicates when the fluid level is low and prompts the user to prepare a new bag of solution.

**Heater On LED:** This is a panel light that indicates when the heater is on. If this light is not on the heater is not maintaining the Pre-set temperature of 100° F.

**Air Connector:** On the Back of panel of the unit a Colder Connector which plugs into a mating connector on the Pressure Enclosure.

**Heater Electrical Connector:** On the back panel of the controller is located a connector used to electrically connect the heater and the temperature sensor to the controller.

### 4.2 DESCRIPTION OF THE OPERATION

Plug the unit into a proper outlet according to the power requirements listed on the label at the back of the unit.

Turn on the power located on the front panel.

Set the Load/Run switch located on the front panel to the LOAD position.

Place a fresh preheated bag of solution into the Solution Carrier and insert this into the pressure enclosure.

The solution bag is connected to the Probe/Tubing set through stick fittings that protrude through an opening at the bottom of the pressure enclosure.

After the stick fitting has been plugged into the bag of solution, lock the piercing in locking collar at the bottom of the pressure enclosure.

Set the LOAD/RUN switch to the run position.

21



#### **4.2 (continued)**

**Set the operating pressure setpoint with the switch located to the right of the LED display located on the front panel.**

**Pushing up on the switch will cause the pressure setting to go up. Pushing down on the switch will cause the setting to down.**

**After a setting has been set for more the 10 seconds the device will store as the last setting and will automatically go to this setting when the is powered up in the next session. Always check the setting before proceeding with procedure.**

**Temperature is automatically set to 100 F°. The heater is designed to maintain the temperature of the irrigation fluid.**

**When the user depresses the fluid control valve on the probe fluid will flow and provide irrigation.**

**When the irrigation fluid level gets low in the bag a FLUID LOW indicator light will come on indicating that level is low in the irrigation bag.**

**When the bag becomes empty set the LOAD/RUN switch to the LOAD position. This will cause the pressure to drop off in the pressure enclosure and a fresh bag of irrigation fluid can be installed.**

**The device will stored the last settings until they are changed by the user.**



## 5.0 USAGE / CARE / WARRANTY

### USAGE RECOMMENDATION

The **THERMO-FLO SYSTEM 3**<sup>TM</sup> was designed for Laparoscopic surgery and is only to be used by licensed, certified surgical personnel.

### RECOMMENDED CARE

Do not leave unit plugged in.

Do not use abrasive cleaning and disinfecting solutions.

### LIMITED WARRANTY

Sanese Medical Corporation ("Sanese") warrants the **THERMO-FLO SYSTEM 3**<sup>TM</sup> to be free from defects in material and workmanship for a period of one (1) year from date of initial purchase. No warranty shall apply if the goods have been damaged by accident, abuse, misuse, or misapplication, or as a result of service or modification by other than a person authorized by Sanese. Sanese's liability under this warranty, and buyer's exclusive remedy is limited to the cost of materials and labor to repair defective goods, or to their replacement, and in no event shall exceed the purchase price. Repair or replacement of defective goods under warranty will be made only upon their return to Sanese after notice to Sanese and Buyer's receipt of shipping instructions. Upon receipt of goods returned under warranty, Sanese will inspect the goods, will notify buyer of the extent of repair or replacement which Sanese will perform under warranty, which shall be conclusive of Sanese's liability.

Sanese is not responsible for incidental consequential damages resulting from the breach of any expressed or implied warranty, including damage to property and to the extent permitted by law, damages for personal injury. The warranty contained herein is in lieu of all warranties, expressed or implied. No statement of any representative shall extend Sanese's liability as herein established or limited.

**SANESE MEDICAL**  
CORPORATION

LOT NO. 0600

**THIS DEVICE IS NON-STERILE AND NEEDS TO BE STERILIZED AND  
CLEANED ACCORDING TO MANUFACTURING SPECIFICATIONS.**

**CONTENTS:** MODEL NO. TFSS - 3 RPR  
1 - Reusable Probe Set

**INDICATION FOR IRRIGATION  
DURING LAPAROSCOPIC PROCEDURES**

**CAUTION:** Federal (U.S.A.) Law restricts this device to sale by,  
or on the order of a physician.  
Recommended for Laparoscopic Irrigation  
in adult and pediatric surgery.

**THIS PRODUCT TO ONLY BE USED FOR  
LAPAROSCOPIC IRRIGATION. DO NOT  
USE FOR HYSTEROSCOPIC DISTENTION.**

*Distributed by:*

**SANESE MEDICAL  
CORPORATION**

490 HORTONWAY ROAD, COLUMBUS, OHIO 43260 (614) 291-8200 FAX (614) 291-1000

LABEL SIZE - 3.25 X 5.50  
UL PDQ2 SPECIFICATIONS  
TEXT COLOR- BLACK  
BACKGROUND COLOR- WHITE

ECN	REV	DESCRIPTION	DATE	DWN		
SANESE MEDICAL		TITLE	TFSS - 3 RPR PRODUCT LABEL			
SHEET 1 OF 1		PART NUMBER	10030	DWN	DATE	
SCALE		FULL	FILE	TFSS3RPR.CDR	MDI	3/1/96

24

**SANESE MEDICAL®**  
CORPORATION

MANUFACTURER OF MEDICAL PRODUCTS

PREMARKET 510(k)  
NOTIFICATION OF CLASS II  
DEVICE. CLASSIFICATION  
REQUESTED BY NOTIFICATION  
85 HET (LAPAROSCOPE  
ACCESSORY); REGULATION NO  
21 C.F.R. § 884.1720

CONFIDENTIAL - 21 CFR § 807.95(b,c)  
CLAIMS MADE  
DO NOT LOG IN PUBLICLY  
AVAILABLE LOG

February 22, 1995


Office of Device Evaluation  
Centers for Devices and Radiological Health  
Document Mail Center (HFZ-41)  
Food and Drug Administration  
9200 Corporate Boulevard  
Rockville, MD 20850

Re: Second Submission of Additional Information in Support  
of Section 510(k) Notification of Sanese Medical  
Corporation's Intention to Market Thermo-Flow System 3  
Laparoscopic Irrigation System and Irrigation  
Aspiration Probe and Tubing  
510(k) No.: K941507/A

Dear Sir/Madam:

Enclosed is an original and two copies of Sanese Medical Corporation's second submission of additional information in support of its Section 510(k) notification of its intention to market its Thermo-Flow System 3 Laparoscopic Irrigation System and Irrigation/Aspiration Probe and Tubing.

Respectfully submitted,



Christopher N. Sanese  
Regulatory Correspondence  
Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, OH 43212  
Telephone: (614) 291-8848  
Fax: (614) 291-1389



**SANESE MEDICAL**<sup>®</sup>  
CORPORATION

MANUFACTURER OF MEDICAL PRODUCTS

PREMARKET 510(k)  
NOTIFICATION OF CLASS II  
DEVICE. CLASSIFICATION  
REQUESTED BY NOTIFICATION  
85 HET (LAPAROSCOPE  
ACCESSORY) REGULATION NO.  
21 C.F.R. § 884.1720

CONFIDENTIAL - 21 CFR § 807.95(b,c)  
CLAIMS MADE  
DO NOT LOG IN PUBLICLY  
AVAILABLE LOG

February 22, 1995

Office of Device Evaluation  
Centers for Devices and Radiological Health  
Document Mail Center (HFZ-41)  
Food and Drug Administration  
9100 Corporate Boulevard  
Rockville, MD 20850

Re: Second Submission of Additional Information in Support  
of Section 510(k) Notification of Sanese Medical  
Corporation's Intention to Market Thermo-Flow System 3  
Laparoscopic Irrigation System and  
Irrigation/Aspiration Probe and Tubing.  
510(k) No.: K941507/A

Dear Sir/Madam:

I.

Pursuant to 21 C.F.R. § 807.87(H) Sanese Medical Corporation (hereinafter "SMC") hereby submits additional information in support of its 510(k) submission under 510(k) No. K941075/A. The additional information provided through this second submission consists simply of promotional material of a predicate device which was inadvertently omitted from the first submission of additional information dated January 13, 1995.

II.

Additional Information

Attached as Exhibit F is a copy of promotional material from Level 1 Technologies Inc. concerning a heated irrigation system indicated for usage during laparoscopic surgery. The attached Exhibit F should be substituted for the present Exhibit F attached to the first submission of additional information in support of the SMC 510(k) notification dated January 13, 1995.

Office of Device Evaluation  
February 22, 1995  
Page - 2 -

The current Exhibit F to the first submission while a similar product also marketed by Level 1 Technologies Inc. is not the same type of heated irrigation system as the SMC device and was erroneously attached to the first submission when it was sent in January 1995.


III.

### CONCLUSION

SMC believes that with the second submission of the additional information described in II above and as supported with the Exhibit F attached hereto, its Thermo-Flow System 3 Laparoscopic Irrigation System and Irrigation/Aspiration Probe and Tubing are substantially equivalent in terms of their safety and effectiveness to the legally marketed predicate devices described in the original 510(k) submission and identified in first and second submissions of additional information. Accordingly, SMC requests that pursuant to 21 C.F.R. § 807.100, FDA issue an order granting SMC premarket clearance to market this device.

Should you have any questions concerning this additional information, please do not hesitate to call or write to me at the number and address set forth below.

Respectfully submitted,



Christopher N. Sanese  
Regulatory Correspondence  
Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, OH 43212  
Telephone: (614) 291-8848  
Fax: (614) 291-1389

**NEW EXHIBIT F**

*28*

  
**LEVEL**

TECHNOLOGIES INC

**Thank you for requesting information on the LEVEL 1 Normothermic Irrigating System. The information has been enclosed.**

**Hypothermia is a major cause of complications. The Level 1 Normothermic Irrigating System provides a high flow supply of body temperature irrigating solution so your patient *stays warm, safe and stable.***


**With Normothermic Irrigating your patients will benefit clinically and your hospital will benefit economically. Normothermic Irrigating makes sense in today's economically sensitive medical environment. A warm stable patient means a fast, uneventful recovery, minimizing valuable nursing time and therefore reducing the overall cost of care.**

**The Level 1 System is easy to use. It requires no change in your present protocol. The automatic air elimination system removes air bubbles, ensuring a clear, unobstructed field of vision. Our new Power Pole lifts up to 12 liters of fluid at the touch of a switch, facilitating leg changes and making flow rate adjustments fast, safe and easy.**

**Let us prove it! The enclosed literature describes the clinical benefits of Normothermic Irrigating in more detail, but you don't have to take our word for it. Try the system on a few cases. The Level 1 Normothermic Irrigating System makes sense, it's practical and your patients will feel the difference!**

**Call us at 1-800-5-LEVEL-1 (1-800-553-8351) to arrange an evaluation.**





# Hypothermia is a dangerous downhill run.

Now, the same system that kept tens of thousands  
of patients warm last year is available for  
*all irrigating applications.*

- **Hysteroscopy**
- **Arthroscopy**
- **Cystoscopy**

In fact, *any* procedure that uses 3 liters  
or more of irrigating fluid.

## The Level 1 Normothermic Irrigating System.

Making surgery safer  
by reducing the incidence  
of hypothermia.




**LEVEL 1**  
TECHNOLOGIES INC. NORMOTHERMIC IRRIGATING

30

# Adaptability.


Need to use a specialized pump for your irrigating procedures? With Level 1, it's easy as 1,2,3.

1. Set up and prime the Level 1 Normo-thermic Irrigating Set.
2. Spike your pump's tubing into the multi-purpose connector.
3. Prime your pump as always.



**The Level 1 Normothermic Irrigating System.  
Adapts to all your needs.**

Questions? Call 1 800 5 LEVEL 1 (1 800 553-8351)



Level 1 Technologies, Inc., 160 Weymouth Street, Rockland, MA 02370  
800 553-8351, 617 878-8011, FAX 617 878-8201

31

Food and Drug Administration  
Center for Devices and  
Radiological Health  
Office of Device Evaluation  
Document Mail Center (HFZ-401)  
9200 Corporate Blvd.  
Rockville, Maryland 20850

January 20, 1995

SANESE MEDICAL CORPORATION  
885 NORTHWEST BLVD.  
COLUMBUS, OH 43212  
ATTN: CHRISTOPHER N. SANESE

510(k) Number: K941507  
Product: THEMO-FLO SYSTEM  
3  
LAPAROSCOPIC/HYS  
TEROSCOPIC

The additional information you have submitted has been received.

We will notify you when the processing of this submission has been completed or if any additional information is required. Please remember that all correspondence concerning your submission MUST be sent to the Document Mail Center (HFZ-401) at the above letterhead address. Correspondence sent to any address other than the one above will not be considered as part of your official premarket notification submission. Because of equipment and personnel limitations we cannot accept telefaxed material as part of your official premarket notification submission, unless specifically requested of you by an FDA official.

The Safe Medical Devices Act of 1990, signed on November 28, states that you may not place this device into commercial distribution until you receive a letter from FDA allowing you to do so. As in the past, we intend to complete our review as quickly as possible. Generally we do so 90 days. However, the complexity of a submission or a requirement for additional information may occasionally cause the review to extend beyond 90 days. Thus, if you have not received a written decision or been contacted within 90 days of our receipt date you may want to check with FDA to determine the status of your submission.

If you have procedural or policy questions, please contact the Division of Small Manufacturers Assistance at (301) 443-6597 or at their toll-free number (800) 638-2041, or contact me at (301) 594-1190.

Sincerely yours,

Marjorie Shulman  
Supervisory Consumer Safety Officer  
Premarket Notification Section  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

32

K941507/S<sup>2</sup>

**SANESE MEDICAL**  
CORPORATION

MANUFACTURER OF MEDICAL PRODUCTS

PREMARKET 510(k)  
NOTIFICATION OF CLASS II  
DEVICE. CLASSIFICATION  
REQUESTED BY NOTIFICATION  
80 HET: REGULATION NO. 21  
C.F.R. § 884.1720

CONFIDENTIAL - 21 CFR § 807.95(b,c)  
CLAIMS MADE  
DO NOT LOG IN PUBLICLY  
AVAILABLE LOG

January 13, 1995

Office of Device Evaluation  
Centers for Devices and Radiological Health  
Document Mail Center (HFZ-41)  
Food and Drug Administration  
9200 Corporate Boulevard  
Rockville, MD 20850

RECEIVED  
17 JAN 20 53  
FDA/CDRH/OGE/DHIC

Re: Submission of Additional Information in Support of  
Section 510(k) Notification of Sanese Medical  
Corporation's Intention to Market Thermo-Flow System 3  
Laparoscopic Irrigation System and Irrigation  
Aspiration Probe and Tubing  
510(k) No.: K941507/A

Dear Sir/Madam:

Enclosed is an original and two copies of Sanese Medical Corporation's submission of additional information in support of its Section 510(k) notification of its intention to market its Thermo-Flow System 3 Laparoscopic Irrigation System and Irrigation/ Aspiration Probe and Tubing.

Respectfully submitted,

Christopher N. Sanese  
Regulatory Correspondence  
Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, OH 43212  
Telephone: (614) 291-8848  
Fax: (614) 291-1389

SS

PREMARKET 510(k)  
NOTIFICATION OF CLASS II  
DEVICE. CLASSIFICATION  
REQUESTED BY NOTIFICATION  
85 HET: REGULATION NO. 21  
C.F.R. § 884.1720

CONFIDENTIAL - 21 CFR § 807.95(b,c)  
CLAIMS MADE  
DO NOT LOG IN PUBLICLY  
AVAILABLE LOG

January 13, 1995

Office of Device Evaluation  
Centers for Devices and Radiological Health  
Document Mail Center (HFZ-41)  
Food and Drug Administration  
9100 Corporate Boulevard  
Rockville, MD 20850

Re: Submission of Additional Information in Support of  
Section 510(k) Notification of Sanese Medical  
Corporation's Intention to Market Thermo-Flow System 3  
Laparoscopic Irrigation System and  
Irrigation/Aspiration Probe and Tubing.  
510(k) No.: K941507/A

Dear Sir/Madam:

I.

Pursuant to the November 17, 1994 letter of Colin M. Pollard, Chief, Obstetrics/Gynecology Devices Branch, Division of Reproductive, Abdominal, Ear, Nose and Throat and Radiological Device, Office of Device Evaluation, Centers for Devices and Radiological Health (copy attached as Exhibit A), Sanese Medical Corporation (hereinafter "SMC") hereby submits pursuant to 21 C.F.R. § 807.87(H) additional information in support of its 510(k) submission under 510(k) No. K941075-A.<sup>1/</sup> The additional information provided through this submission consists of certain tests, some of which were in existence at the time of the original submission but were deemed unnecessary for inclusion in the notification, an amended summary of safety and effectiveness information, software validation results and modified labeling. To aid in understanding the submission, the requests for

<sup>1/</sup> Pursuant to SMC's written request, it was granted until February 17, 1995 to submit additional information in response to the request. (See Exhibit B.)

Office of Device Evaluation  
January 13, 1995  
Page - 2 -

additional information set forth in Mr. Pollard's November 17, 1994 letter are set out in their entirety with SMC's responses following each numbered request and/or its subpart.

II.

b(4)CCI



Office of Device Evaluation  
January 13, 1995  
Page - 3 -

b(4)CCI



Office of Device Evaluation  
January 13, 1995  
Page - 4 -

b(4)CCI





Office of Device Evaluation  
January 13, 1995  
Page - 5 -

b(4)CCI



Office of Device Evaluation  
January 13, 1995  
Page - 6 -

b(4)CCI



39

b(4)CCI



III.

b(4)CCI



40

Office of Device Evaluation  
January 13, 1995  
Page - 8 -

b(4)CCI



A

42



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
1380 Piccard Drive  
Rockville MD 20850

Mr. Christopher N. Sanese  
Regulatory Correspondent  
Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, Ohio 43212

NOV 17 1994

Re: K941507/A  
Thermo-Flo System 3  
Dated: July 5, 1994  
Received: July 6, 1994

Dear Mr. Sanese:

We have reviewed your Section 510(k) notification of intent to market the device referenced above. We cannot determine if the device is substantially equivalent to a device marketed prior to May 28, 1976, the enactment date of the Medical Device Amendments, based solely on the information you provided.

(b)(4) Trade Secret Process - Product Specification

(b)(4) Trade Secret Process - Product Specification

UB

(b)(4) Trade Secret Process - Product Specification



614

We believe that this information is necessary for us to determine whether or not this device is substantially equivalent to a legally marketed predicate device with regard to its safety and effectiveness.

You may not market this device until you have provided adequate information described above and required by 21 CFR 807.87(f) and (h), and you have received a letter from FDA allowing you to do so. If you market the device without conforming to these requirements, you will be in violation of the Federal Food, Drug, and Cosmetic Act (Act). You may, however, distribute this device for investigational purposes to obtain clinical data if needed to establish substantial equivalence. Clinical investigations of this device must be conducted in accordance with the investigational device exemption (IDE) regulations.

The requested information should reference your above 510(k) number and should be submitted in duplicate to:

Food and Drug Administration  
Center for Devices and  
Radiological Health  
Document Mail Center (HFZ-401)  
9200 Corporate Boulevard  
Rockville, Maryland 20850

If the information is not received within 30 days, we will consider your premarket notification to be withdrawn and your submission will be deleted from our system. If you submit the requested information after 30 days it will be considered and processed as a new 510(k); therefore, all information previously submitted must be resubmitted so that your new 510(k) is complete.

If you have any questions concerning the contents of this letter, please contact Donna-Bea Tillman at (301) 594-1180. If you need information or assistance concerning the IDE regulations, please contact the Division of Small Manufacturers Assistance at their toll free number (800) 638-2041 or at (301) 443-6597.

Sincerely yours,

*Colin M. Pollard*

Colin M. Pollard  
Chief, Obstetrics/Gynecology Devices Branch  
Division of Reproductive,  
Abdominal, Ear, Nose and Throat,  
and Radiological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

US



Page 4 - Mr. Christopher N. Sanese

bcc: HFZ-401 DMC  
HFZ-404 510(k) Staff  
HFZ- Division  
D.O.

DBTillman  
C:\WP51\FILES\K941507.AI  
Draft 11-15-94 ens

46

B

47

12-9-94

Colin M. Pallard  
Food and Drug Administration  
Center for Devices and  
Radiological Health  
Document Mail Center (HFZ-401)  
9200 Corporate Boulevard  
Rockville, Maryland 20850

RE:K941507 / A  
Thermo - Flo System 3  
Dated: July 5, 1994  
Received: July 6

Dear Mr. Pallard,

This letter is to formally request a 60 day extension of our December 17, delivery compliance as requested of your letter of November 17, 1994. This request was made verbally to Donna -Bea Tillman in our phone conversation of December 5, 1994.

Thank you very much for your consideration of this request.

Sincerely,

Christopher N. Sanese  
Regulatory Correspondent  
SANESE MEDICAL CORPORATION

ENC: 4 pages

CNS/III

CC: John Jevicky, Esq.  
DINSMORE & SHOHL

48

Food and Drug Administration  
Center for Devices and  
Radiological Health  
Office of Device Evaluation  
Document Mail Center (HFE-401)  
9200 Corporate Blvd.  
Rockville, Maryland 20850

December 15, 1994

SANESE MEDICAL CORPORATION  
885 NORTHWEST BLVD.  
COLUMBUS, OH 43212  
ATTN: CHRISTOPHER N. SANESE

510(k) Number: K941507  
Product: THEMO-FLO SYSTEM  
3  
LAPAROSCOPIC/HYS  
TEROSCOPIC

Extended Until: 17-FEB-95

Based on your recent request, an extension of time has been granted for you to submit the additional information we requested.

If the additional information is not received by the "Extended Until" date shown above your premarket notification will be considered withdrawn.

If you have procedural or policy questions, please contact the Division of Small Manufacturers Assistance at (301) 443-6597 or at their toll-free number (800) 638-2041, or contact me at (301) 594-1190.

Sincerely yours,

*Marjorie Shulman*

Marjorie Shulman  
Supervisory Consumer Safety Officer  
Premarket Notification Section  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

49

C

SD

AMENDED SECTION 510(K) SUMMARY OF SAFETY AND  
EFFECTIVENESS INFORMATION FOR SMC'S  
THERMO-FLOW SYSTEM 3 LAPAROSCOPIC  
IRRIGATION SYSTEM AND  
DISPOSABLE IRRIGATION/ASPIRATION  
PROBE AND TUBING

Pursuant to Section 12(a) of the Safe Medical Devices Act of 1990, 21 U.S.C. § 360c(i)(3)A and 21 C.F.R. § 807.92, the following amended summary of information concerning the safety and effectiveness of SMC's Thermo-Flow System 3 Laparoscopic Irrigation System and Disposable Irrigation/Aspiration Probe and Tubing is hereby submitted.

1. Name and address of device manufacturer submitting 510(k) notification:

Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, OH 43212  
(614) 291-8848

2. Regulatory correspondent.

Christopher N. Sanese  
Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, OH 43212  
(614) 291-8848

3. Date amended summary was prepared: January 10, 1995

4. Names of Devices:

- (i) Proprietary Names:

- (a) Thermo-Flow System 3 Laparoscopic Irrigation System\*
- (b) Sanese Disposable Irrigation/Aspiration Probe and Tubing™

- (ii) Common, usual names:

- (a) Laparoscopic irrigation pump system
- (b) Disposable aspiration/irrigation probe and tubing.

- (iii) Classification Names: Gynecologic Laparoscope and Accessories, 80 HET: 21 C.F.R. § 884.1720.

5. Predicate devices to which SMC is claiming substantial equivalence:

51

SMC claims substantial equivalence for purposes of section 510(k) only to:

A. Irrigation Pump System:

- (1) The Niagra™ High-Flow Irrigator System, 510(k) Number K924530 by Cabot Medical, 2021 Cabot Boulevard, West Langhorn, PA 19047 (800-523-6078).
- (2) Level One Technologies, Inc., Labroscopic Irrigation System, 160 Weymouth Street, Rockland, MA 02370 (800-553-8351).

B. Disposable Probe and Tubing:

(1) Probe:

Corson™ Disposable Suction Irrigation Probe by Cabot Medical, 2021 Cabot Boulevard, West Langhorn, PA 19047 (800-523-6078).

(2) Tubing:

Core Dynamics Inc. Disposable Irrigation/Suction Tubing Set (510(k) Number K910246) by Core Dynamics, Inc., 9951 Atlantic Boulevard, Jacksonville, Florida 3225 (904-727-0910).

6. Description of the devices and primary principles of operation:

A. Description of irrigation pump system.

This device consists of a fabric enclosure (pressure envelope) in which a bag of irrigation fluid is placed. The fabric enclosure contains an inflatable bladder that is pneumatically inflated to pressurize the fluid in the irrigation fluid bag. Inflation of the bladder is provided by an air compressor which has a maximum voltage leakage current of 100 microamps. This complies with Underwriters Laboratory Standard 544.

The pressure envelope also contains a heating element to assist in maintaining the temperature of 100° F. of preheated irrigation solution while being used in the pressure envelope. There is also a thermal regulator to control the temperature of the heating element. The heating element and the thermal regulator are powered by an Underwriter Laboratory approved 24 volt electrical power supply.

The air compressor and heating element are controlled by a microcontroller. These components are controlled by the micro-controller through activating solid state switches which control electrical power to the pump, pump solenoids and the heater.

The pressurized preheated fluid is then introduced into the laparoscopic surgery site through a medical grade plastic tube and probe that are discussed below.

B. Description of disposable probe and tubing:

The probe consists of a medical grade plastic, activation valve and metal shaft that will deliver and aspirate fluids to and from an operative site.

The tubing consists of two 15 foot length medical grade surgical tubing that are approximately 10 millimeters in diameter which would be connected to the probe and the irrigation solution bag and the pump when used for irrigation or to the probe and the suction unit when used for aspiration. The probe and tubing are made of medical grade plastic and metal that have been determined by FDA to be safe and effective for this intended usage.

C. Primary Principles of Operation:

(i) Irrigation Pump System:

The primary principles of operation are pressurization through an electrically powered pneumatic compressor for pressurized irrigation fluid flow and thermal transfer through an electrically-powered thermo-conducted heat transfer unit to maintain preheated fluid.

(ii) Probe and Tubing:

With respect to the aspiration/irrigation use of the probe and tubing, the principles of operation are vacuum and pressurization.

7. Statement of the purpose for which the device will be recommended:

These devices are intended to be used accessories to laparoscopic surgical equipment by providing irrigation and aspiration during adult and pediatric laparoscopic surgery.

8. Biocompatibility and effectiveness of the device:



(i) Biocompatibility:

The only aspect of the device that comes in contact with the patient is the distal end of the probe. The probe and the probe tubing have been subject to biocompatibility tests and are safe and effective for such usage. The irrigation pump has a maximum voltage leakage current of 100 microamps in accordance with U.L. Standard 544. Performance tests of the pump demonstrate that it is safe and effective for its intended use. Also, all components of the pump meet U.L. Standards.

Respectfully submitted,

Christopher N. Sanese  
Regulatory Correspondent  
Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, OH 43212  
(614) 291-8848

D

SS

**Preliminary Discussion:**

b(4)Trade Secret Process-Product Specs



**Conditions:**

b(4)Trade Secret Process-Product Specs



**Data:**

b(4)Trade Secret Process-Product Specs



SP

**Results:**

b(4)Trade Secret Process-Product Specs



57

b(4)Trade Secret Process-Product Specs

A large black rectangular redaction box covering the top portion of the page content.

#### **4 Numerical Specifications**

b(4)Trade Secret Process-Product Specs

A large black rectangular redaction box covering the middle portion of the page content.

#### **5 User Interface**

##### **5.1 Description of Controls**

b(4)Trade Secret Process-Product Specs

A large black rectangular redaction box covering the bottom portion of the page content.

58

E

59

Sanese Medical Test Report:  
Test Number: eo001-12 REV A  
Subject : Testing of Machine Setting Accuracy  
Written: Marc D. Taylor  
Date: 12/12/94  
Page 1 of 2

*W/T D 7/92 12/12/94*

**Objective:**

b(4)Trade Secret Process-  
P d tS

**Conditions:**

b(4)Trade Secret Process-Product Specs

**Results:**

b(4)Trade Secret Process-Product Specs

**Discussion:**

b(4)Trade Secret Process-Product Specs

*60*

Sanese Medical Test Report:  
Test Number: eo001-12 rev a  
Subject: Testing of Machine Setting Accuracy  
Written: Marc D. Taylor  
Date: 12/12/94  
Page 2 of 2

*Marc D. Taylor 12/12/94*

## GRAPH OF AIR AND FLUID PRESSURES

b(4)Trade Secret Process-Product Specs



*61*



Sanese Medical Test Report:  
Test Number: eo001-16 REV A  
Subject : Flow Rate Test  
Written: Marc D. Taylor  
Date: 12/12/94  
Page 1 of 2

**Objective:**

[Redacted] b(4)Trade Secret [Redacted]

**Conditions:**

b(4)Trade Secret Process-Product Specs  
[Redacted]

**Test:**

b(4)Trade Secret Process-Product Specs  
[Redacted]

**Results:**

b(4)Trade Secret Process-Product Specs  
[Redacted]

**Discussion:**

b(4)Trade Secret Process-Product Specs  
[Redacted]

602

**Sanese Medical Test Report**  
**Test Report Number: eo001-16 revA**  
**Written By: Marc D. Taylor**  
**Date: 12/12/94**  
**Page 2 of 2**

b(4)Trade Secret Process-Product Specs




63

F

64

K911383 - DGRD

# Hypothermia makes TURP a dangerous downhill run.



Ask your Recovery Room staff.  
TURP patients are cold!

Even *lukewarm* irrigating during TURP makes core temperature loss unavoidable. This is significant because hypothermia increases cardiac stress, elevates bleeding times and lengthens the recovery period.

### Eliminate the risk of hypothermia.

Normothermic Irrigating lowers the risk of serious hypothermic complications, reduces bleeding and maintains stable cardiac function. With the Level 1 Normothermic Irrigating System, regardless of the length of the procedure or the volume of irrigant used, you can be confident your patient will leave your O.R. warm. *In fact, the more fluid used, the warmer your patient will stay!*

### Specifically for Urology.

The H-600 Normothermic Irrigating System provides a continuous, high flow supply of body temperature irrigating fluid so your patient *stays* warm, safe and stable. An in-line Air Eliminator assures a clear, bubble-free field of vision. And, our exclusive Power Pole™ lifts up to four three-liter bags at the touch of a switch... bag changes and height adjustments become quick, safe and easy.

Evaluate the H-600 Normothermic Irrigating System. Your patients will feel the difference!

An Eliminator  
insure's bubble  
free field of  
vision.



Level 1 Technologies, Inc., 160 Weymouth Street, Rockland, MA 02370, 1 (800) 5 LEVEL 1, FAX 617 878 8201  
©1992 Level 1 Technologies, Inc.

Handwritten initials 'US'.

K861477 - DCRND / K873435 - GU



**LEVEL**

TECHNOLOGIES INC.

## **NORMOTHERMIC IRRIGATING DURING TURP?**

**-Room temperature and lukewarm irrigating solutions cause a serious loss of core temperature.<sup>1,2,3,4</sup>**

**-This temperature loss causes significant cardiovascular stress<sup>1,2,7</sup> and increases the risk of additional hypothermic complications.<sup>1,2,3,4,5,6,7</sup>**

**-Normothermic irrigating maintains core temperature, reduces bleeding and maintains safe, stable cardiac function.<sup>2,5,6</sup>**

**The Level 1 Normothermic Irrigating System eliminates hypothermic complications and provides for "a less stressful and safer operation" <sup>5</sup>**

## **NORMOTHERMIC IRRIGATING DURING TURP!**

Call 800-5-LEVEL-1 to arrange an evaluation at your facility.

- 1. Haemodynamic evidence for cardiac stress during transurethral prostatectomy, J.W.H. Evans, et al, *British Medical Journal* 1992;304:666-671.**
- 2. Haemodynamic and cardiological responses to TURP, standard vs. isothermic techniques, J.W.H. Evans, et al, *presented to the British Society of Urological Surgeons*, June, 1991**
- 3. Hypothermia during transurethral resection of the prostate, Carpenter, *Urology*, 2/84**
- 4. Body temperature changes during prostatic resection as related to the temperature of the irrigating solution, Allen, *The Journal of Urology*, 10/73**
- 5. The effect of warm irrigation on blood loss during transurethral prostatectomy under spinal anesthesia, Heathcote and Dyer, *British Journal of Urology*, 4/86**
- 6. The effects of warm irrigating fluid during and after transurethral prostatectomy, Ogura, et al, *Clinical Therapeutics*/Vol 10, special issue, 1988**
- 7. Hypothermia is An Independent Predictor Of Postoperative Myocardial Ischemia, S.M. Frank, et al, *Anesthesiology*, V77, No 3A, September, 1992**

*\*Summaries of the above studies are attached. The complete articles are available in your library or on request from*

**LEVEL 1 TECHNOLOGIES**

**Haemodynamic evidence for cardiac stress during transurethral prostatectomy,**  
J. W. H. Evans, et al, *British Medical Journal* 1992;304: 666-671.

This paper examines changes in cardiac output, stroke volume, systemic vascular resistance, and core temperature between a group of patients undergoing routine TURP (intermittent flow technique) and a control group, undergoing either hernia repair or testicular exploration. The same anesthetic protocol was used in both groups. Both groups compare well in terms of age of patient and length of procedure.

**RESULTS**

**CONTROL GROUP**

- No significant changes in hemodynamic parameters measured
- Core temperature was stable

**ROUTINE TURP GROUP**

- Systemic vascular resistance rose significantly from baseline compared to the control group (+49% vs Control -10%)
- Cardiac output fell significantly (-21% vs Control +12%)
- Stroke volume also fell significantly (-16% vs Control +26%)
- Core temperature fell significantly compared to the control group

**AUTHOR'S CONCLUSIONS**

**"Important haemodynamic(sic) disturbances were identified during routine apparently uneventful transurethral prostatectomy but not during control procedures."**

**"The haemodynamic differences between groups must be the result of a factor peculiar to transurethral prostatectomy which was present early in the procedure."**

**"The most remarkable difference between the two groups was the rapid central cooling during transurethral prostatectomy"**

**Haemodynamic and Cardiological Responses to TURP: Standard vs. Isothermic Techniques, J.W.H. Evans et al, presented to the British Society of Urological Surgeons, June, 1991.**

This study compares cardiovascular and core temperature changes in patients undergoing TURP using room temperature irrigant, TURP using normothermic (body temperature) irrigant, and a control group undergoing non-endoscopic surgery. The same anesthetic protocol was used in all groups. The duration of surgery was comparable in all groups.

**RESULTS**

**Non-endoscopic surgery (no bladder irrigation)**

- No significant changes in heart rate, stroke volume, cardiac output, or vascular resistance
- Core temperature was stable

**TURP with normothermic irrigant**

- No significant changes in the above hemodynamic parameters or core temperature

**TURP with room temperature irrigant**

- Mean arterial pressure and vascular resistance rose significantly (over 55%)
- Stroke volume and cardiac output fell significantly (-20% and -25% respectively)
- Core Temperature dropped significantly
- 20% of the room temperature TURP group became hypothermic during surgery (core temperature < 35°C). [Previous study documents that core temperature continues to drop in recovery...Allen, *The Journal of Urology*, 10/73]

**AUTHORS' CONCLUSION**

- "The hemodynamic responses during standard TURP [room temperature irrigant] were the result of falls in core temperature..."
- "Iatrogenic falls in core temperature represent a cardiovascular stress to all TURP patients. Furthermore, low core temperatures have many other detrimental physiological and metabolic effects, especially in the elderly."
- "The prevention of falls in core temperature during TURP is a simple, practical, and effective preventative measure."

108

**Hypothermia During Transurethral Resection of Prostate, Carpenter, Urology 2/84**

*Author reviews previous articles and conducts his own research on using room temperatures vs body temperature irrigating fluids.*

**Previous Study**

Room temperature solutions drop body temperature 2oF per hour.

120 min. latter body temperature still dropping, the study was concluded at this time.

Warmed solutions: no drop in temperature; no increase in bleeding.

**Carpenter's Results**

Average operating time 90 minutes

Average patients temperature dropped for 180 minutes.

Temperatures returned to normal 7 hours after onset of anesthesia.

Several patients remained cold for up to 15 hours.

**Comments**

Elderly may have defective thermoregulation. This can result in a lack of shivering to protect core temperatures.

Cardiovascular system especially sensitive to hypothermia.

Hypothermic heart relatively unresponsive to atropine, electropacing or counter shock.

Turp with room temperature relatively safe procedure.

Geriatric group may have pre-existing cardiovascular disease thus possibility of complications is always present.

**Authors Conclusions**

*"The concept that mild hypothermia is beneficial in regard to intraoperative blood loss would seem to be greatly outweighed by the potentially serious cardiovascular complications which might occur as a result of hypothermia."*

69



**BODY TEMPERATURE CHANGES DURING PROSTATIC RESECTION AS RELATED TO THE TEMPERATURE OF THE IRRIGATING SOLUTION**, Allen, The Journal of Urology, 10/73

*This study looks at the effect on body temperature during TURP using irrigating solutions at three different temperatures;*

**1. Cool (68-84°F) - 2. Lukewarm (85-94°F) - 3. Warm (95-100°F)**

**Author's Observations**

- Due to conventional O.R. temperatures patients tend to have a decline in body temperature even before resection has begun.
- This decline accelerates upon the institution of cool OR lukewarm irrigation of the bladder.
- The cooling effect was only *SLIGHTLY* less when the irrigant was lukewarm. There was an "invariable fall in body temperature when cool OR lukewarm irrigation solutions were used." (emphasis ours)
- In the warm series, no drop in body temperature occurred.
- A fall in body temperature of as little as 2/3°C may cause shivering.
- In the cardiac patient with ischemic heart disease, shivering may be sufficient to induce arrhythmias, angina or even myocardial infarction.

**Author's Conclusions**

*"The use of warm irrigating solutions seemed to have considerable superiority over other methods of maintaining body temperature."*

*Use of warm solutions makes the patient and surgeon more comfortable.*

*No adverse effects were seen with the use of warm solutions, but cool solutions caused shivering in 16 per cent of the patients."*

70

**The effect of Warm Irrigation on Blood Loss during Transurethral Prostatectomy under Spinal Anesthesia, Heathcote and Dyer, British Journal of Urology, 4/86**

*This study describes the effect of warm irrigating fluid on blood loss during TURP.*

**Authors' comments:**

- Hypothermia is a common problem during TURP.
- Hypothermia causes a risk of cardiovascular complications.
- Cold irrigating fluids cause discomfort for patient and operator.
- With cooled irrigating fluids, shivering and exaggerated temperature decrements up to 6°C have been reported.
- Warm irrigating fluids are more comfortable for the operator and the increased pliability of the tubing made the resectoscope easier to handle.

**Authors' Conclusion:**

*"By using warmed irrigating fluids during TURP one can decrease heat loss and shivering in the patients. Warmed irrigating fluids do not increase blood loss during TURP and their use leads to a less stressful and safer operation."*

71

**The Effects of Warm Irrigating Fluid during and after Transurethral Prostatectomy, Ogura et al, Clinical Therapeutics/Vol 10, special issue, 1988**

*This study compares warm (37°C) versus room temperature irrigation in TURP patients to determine effects on total blood loss.*

**Previous literature documents:**

- Irrigation during TURP with room temperature fluids causes a decrease in core temperature.
- Intraoperative hypothermia can seriously affect the cardiovascular system in elderly patients.

**This study shows:**

- Normothermic Irrigating fluid group showed significantly less blood loss than the control group.
- Far less shivering during and after surgery was noted in the warmed-fluid group (8% vs. 35% in the room temperature group).

**Authors' conclusions:**

*"The results of the present investigation demonstrate that transurethral resection is safer when irrigating fluid is warmed to body temperature (37°C). Warmed irrigating fluid appeared to reduce blood loss locally and prevented hypothermia and associated complications systemically."*

72

# HYPOTHERMIA IS AN INDEPENDENT PREDICTOR OF POSTOPERATIVE MYOCARDIAL ISCHEMIA

S.M. Frank, et al, *Anesthesiology*, V77, No 3A, September, 1992

It has been previously suggested that shivering may be detrimental in patients with coronary artery disease, but no documentation exists to support this hypothesis. This abstract is from a study conducted at Johns Hopkins Hospital. It examines the relationship between inadvertent hypothermia and cardiac morbidity.

## METHODS

Patients' temperatures were recorded upon admission to the ICU (hypothermia  $<35^{\circ}\text{C}$ , normothermia  $\geq 35^{\circ}\text{C}$ ), shivering was determined by physical examination and ischemia was identified via 24 hour holter monitoring.

## RESULTS

The incidence of myocardial ischemia was almost **THREE TIMES GREATER** in the hypothermic group (36% vs 13% for the normothermic group)

Shivering was *not* a predictor of ischemia, "the only independent predictor of ischemia was temperature on arrival in the ICU."

The odds of developing myocardial ischemia increased by 82% for each degree C decrease in core temperature.

## AUTHORS' CONCLUSIONS

"The incidence of postoperative myocardial ischemia was significantly greater in the hypothermic patients"

"Contrary to expectation, the incidence of ischemia was not related to shivering."

"... patients at risk for perioperative myocardial ischemia should receive special attention regarding the maintenance of body temperature."

G

74

**Objective:**

b(4)Trade Secret Process-Product Specs



**Conditions:**

b(4)Trade Secret Process-Product Specs



**Results:**

b(4)Trade Secret Process-Product Specs



**Discussion:**

b(4)Trade Secret Process-Product Specs



75

**Sanese Medical Test Report:**

**Test Number: 0001-06 REV A**

**Subject : Production - Line Dielectric Voltage - Withstand Test  
Production - Line Grounding - Continuity Test**

**Written: Marc D. Taylor**

**Date: 12/12/94**

**Page 1 of 1**

*Marc D. Taylor 12/12/94*

**Objective:**

b(4)Trade Secret Process-Product Specs

**Conditions:**

b(4)Trade Secret Process-Product Specs

**UL Code Table Reference:**

b(4)Trade Secret Process-Product Specs

**Results:**

b(4)Trade Secret Process-Product Specs

**Discussion:**

b(4)Trade Secret Process-Product Specs

*76*

Sanese Medical Test Report:  
Test Number: eo001-13 REV A  
Subject : 35 LB Pull Test of Cord and Strain Relief (7.2.3.4 of the UL 544 Code)

Written: Marc D. Taylor  
Date: 12/12/94  
Page 1 of 1

*Marc D Taylor 12/12/94*

**Objective:**

[Redacted] b(4)Trade  
S t

**Conditions:**

b(4)Trade Secret Process-Product Specs  
[Redacted]

**Test:**

b(4)Trade Secret Process-Product Specs  
[Redacted]

**Results:**

b(4)Trade Secret Process-Product Specs  
[Redacted]

**Discussion:**

b(4)Trade Secret Process-Product Specs  
[Redacted]

*77*



H

78

Descriptions of Functions: (Structure Chart on page 3)

b(4)Trade Secret Process-Product Specs



79

b(4)Trade Secret Process-Product Specs



80

b(4)Trade Secret Process-Product Specs



81

II

82

0 Sub Documents

b(4)Tr  
d

1 Purpose

b(4)Trade Secret Process-Product Specs

2 System Overview

b(4)Trade Secret Process-Product Specs

3 Hardware Overview

b(4)Trade Secret Process-Product Specs

83

## 4 Software Requirements

### 4.1 Software Requirements for the User Interface

b(4)Trade Secret Process-Product Specs



### 4.2 Software Requirements for the Hardware Interface

b(4)Trade Secret Process-Product Specs



84

b(4)Trade Secret Process-Product Specs



85



b(4)Trade Secret Process-Product Specs



86

b(4)Trade Secret Process-Product Specs



**5 Software Specifications:**


b(4)Trade Secret Process-Product Specs



b(4)Trade Secret Process-Product Specs



88



89

b(4)Trade Secret Process-Product Specs



90

b(4)Trade Secret Process-Product Specs



91

b(4)Trade Secret Process-Product Specs



92

b(4)Trade Secret Process-Product Specs



93



K

94

**Sanese Medical Controlled Document**  
**Software Engineering Routine**  
**Written: K.Kelly**  
**Last revised on ECN: N/A**

**No: er004-00 Rev: A**  
**Page 1 of 4**  
**Released: 11/6/94**

---

b(4)Trade Secret Process-Product Specs



as

**Sanese Medical Controlled Document**  
**Software Engineering Routine**  
**Written: K.Kelly**  
**Last revised on ECN: N/A**

**No: er004-00 Rev: A**  
**Page 2 of 4**  
**Released: 11/6/94**

(b)(4)Trade Secret Process-Product Specs



*Handwritten signature or initials*

**Sanese Medical Controlled Document**  
**Software Engineering Routine**  
**Written: K.Kelly**  
**Last revised on ECN: N/A**

**No: er004-00 Rev: A**  
**Page 3 of 4**  
**Released: 11/6/94**

---

b(4)Trade Secret Process-Product Specs



*CP 1*

**Sanese Medical Controlled Document**  
**Software Engineering Routine**  
**Written: K.Kelly**  
**Last revised on ECN: N/A**

**No: er004-00 Rev: A**  
**Page 4 of 4**  
**Released: 11/6/94**

---

b(4)Trade Secret Process-Product Specs



98

L

99

b(4)Trade Secret Process-Product Specs



100

b(4)Trade Secret Process-Product Specs



101



b(4)Trade Secret Process-Product Specs



102

b(4)Trade Secret Process-Product Specs



103

b(4)Trade Secret Process-Product Specs



104

b(4)Trade Secret Process-Product Specs



105

b(4)Trade Secret Process-Product Specs



106

b(4)Trade Secret Process-Product Specs



107

b(4)Trade Secret Process-Product Specs



108

b(4)Trade Secret Process-Product Specs



104



b(4)Trade Secret Process-Product Specs



110

M

|||

b(4)Trade Secret Process-Product Specs



112

2

113

# **SANESE MEDICAL CORPORATION, INC.**

## **REUSABLE PROBE STERILIZATION PROCEDURE**

### Preparation for cleaning

1- Dismantle.

- a) Remove housing screws.
- b) Remove valve gates.
- c) Remove springs.
- d) Remove seals.

2- Immerse housing screws, valve gates, springs and, seals in instrument disinfectant.

3- Wipe down the exterior with a lint free cloth soaked in instrument disinfectant.

4- Rinse interior of the probe by introducing the same disinfectant into the embodiment of the probe and allowing the disinfecting solution to run through the irrigation tube.

### Manual cleaning

1- Completely immerse the dismantled probe into the disinfecting solution. Ensure that all of the air is forced out of the probe by tilting it slightly and moving it to and fro in the solution until all air bubbles have been removed.

2- Do not use abrasives for cleaning the probe. Use only swabs, soft bristle brushes or, cleaning guns on the probe.

3- Check the instrument for adhering protein prior to sterilization. If remnants of protein are discovered resoak in disinfecting solution and repeat step #4.

774

Page 2  
Probe Sterilization

Ultrasonic cleaning

- 1- This probe may be cleaned ultrasonically.
- 2- Check the instrument for adhering protein prior to sterilization. If remnants of protein are discovered reclean ultrasonically or soak the probe in a disinfecting solution and wipe with a swab or a soft bristle brush until the probe is free of protein and, repeat ultrasonic cleaning.


Autoclave sterilization

- 1- Autoclave using hospital standardized procedures. Normally autoclaving is performed with saturated steam at 134 degrees C.
- 2- Steam used for sterilization must be free from any contamination and should neither impede the process nor damage the sterilizer or the probe. To guarantee this meet the recommendations of pr EN 285 regarding the quality of water in the tank as well as the condensate.
- 3- After autoclave sterilization the probes must be stored dry until ready for reuse.

Reassemble and test for leaks

- a) Check seals for wear. If worn replace with new seals provided by the manufacturer. If seals show no signs of wear, proceed to "b."
- b) Replace seals into probe housing.
- c) Replace springs.
- d) Replace valve gates.
- e) Replace housing screws.

115



116



LOT NO. 0333

CONTENTS:

1 -MODEL TFS3 LAPAROSCOPIC IRRIGATION SYSTEM

CAUTION:

Federal (U.S.A.) Law restricts this device to sale by,  
or on the order of a physician.

Non - Sterile : Recommended for Laparoscopic Irrigation  
in adult and pediatric surgery

**SANESE MEDICAL**  
CORPORATION

*Distributed by:*

600 NORTH AVENUE BLVD. COLLAMER, OHIO 43022 (614) 291-2200 FAX (614) 291-2200

LABEL SIZE - 3.25 X 5.50  
UL PDQ2 SPECIFICATIONS  
TEXT COLOR- BLACK  
BACKGROUND COLOR- WHITE

ECN	REV	DESCRIPTION	DATE	DWN	
SANESE MEDICAL		TITLE	TFS3 PRODUCT LABEL		
		PART NUMBER	10020	DWN	DATE
SHEET 1 OF 1	SCALE	FULL	FILE	TFS3.CDR	MDI 12/12/94

117



LOT NO. 0500  
 FOR SINGLE USE ONLY CONTENTS  
 IN UNOPEN, UNMAGED PACKAGES

**CONTENTS:** MODEL NO. TFS3 - 3 WP

1 - Single use, disposable 15' irrigation / aspiration tube set with probe  
 INDICATION FOR IRRIGATION  
 DURING LAPAROSCOPIC PROCEDURES

**CAUTION:** Federal (U.S.A.) Law restricts this device to sale by,  
 or on the order of a physician.  
 Recommended for Laparoscopic irrigation  
 In adult and pediatric surgery.

THIS PRODUCT TO ONLY BE USED FOR  
 LAPAROSCOPIC IRRIGATION. DO NOT  
 USE FOR HYSTEROSCOPIC DISTENTION.

Distributed by:

**SANESE MEDICAL**  
 CORPORATION

200 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-2242 FAX (614) 291-2222

LABEL SIZE - 3.25 X 5.50  
 UL PD62 SPECIFICATIONS  
 TEXT COLOR- BLACK  
 BACKGROUND COLOR- WHITE

ECN	REV	DESCRIPTION	DATE	DWN
SANESE MEDICAL		TITLE	TFS3 - 3 WP PRODUCT LABEL	
		PART NUMBER	10021	DWN
		DATE		
SHEET 1 OF 1	SCALE	FULL	FILE	TFS33WP_CDR MDI 12/12/94

118

LOT NO. 0500  
 FOR SINGLE USE ONLY CONTENTS  
 IN UNOPEN, UNBAMGED PACKAGES

**CONTENTS:** MODEL NO. TFS3 - 3 WOP  
 1 - Single use, disposable 15' Irrigation / aspiration tube set  
 INDICATION FOR IRRIGATION  
 DURING LAPAROSCOPIC PROCEDURES

**CAUTION:** Federal (U.S.A.) Law restricts this device to sale by,  
 or on the order of a physician.  
 Recommended for Laparoscopic Irrigation  
 In adult and pediatric surgery.

THIS PRODUCT TO ONLY BE USED FOR  
 LAPAROSCOPIC IRRIGATION. DO NOT  
 USE FOR HYSTEROSCOPIC DISTENTION.

*Distributed by:*

**SANESE MEDICAL**  
 CORPORATION

888 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 391-6200 FAX (614) 391-6200

LABEL SIZE - 3.25 X 5.50  
 UL PDQ2 SPECIFICATIONS  
 TEXT COLOR- BLACK  
 BACKGROUND COLOR- WHITE

ECN	REV	DESCRIPTION	DATE	DWN	
SANESE MEDICAL			TITLE TFS3 - 3 WOP PRODUCT LABEL		
		PART NUMBER 10022	DWN	DATE	
SHEET 1 OF 1	SCALE	FULL	FILE	TFS33WOP.CDR	MDI 12/12/94

119

LOT NO. 0500  
 FOR SINGLE USE ONLY CONTENTS  
 IN UNOPEN, UNBAMGED PACKAGES

**CONTENTS:** MODEL NO. TFS3 - N- 3 WOP  
 10 - Single use, disposable 15' irrigation / aspiration tube set  
 INDICATION FOR IRRIGATION  
 DURING LAPAROSCOPIC PROCEDURES

**CAUTION:** Federal (U.S.A.) Law restricts this device to sale by,  
 or on the order of a physician.  
 Recommended for Laparoscopic Irrigation  
 In adult and pediatric surgery.

THIS PRODUCT TO ONLY BE USED FOR  
 LAPAROSCOPIC IRRIGATION. DO NOT  
 USE FOR HYSTEROSCOPIC DISTENTION.

Distributed by:

**SANESE MEDICAL**  
 CORPORATION

888 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-2040 FAX (614) 291-4888

LABEL SIZE - 3.25 X 5.50  
 UL PDQ2 SPECIFICATIONS  
 TEXT COLOR- BLACK  
 BACKGROUND COLOR- WHITE

ECN	REV	DESCRIPTION	DATE	DWN
SANESE MEDICAL		TITLE	TFS3 - N- 3 WOP CASE LABEL	
		PART NUMBER	10023	DWN DATE
SHEET 1 OF 1	SCALE	FULL	FILE	TFS3NWOP.CDR MDI 12/12/94

120

LOT NO. 0500  
 FOR SINGLE USE ONLY CONTENTS  
 IN UNOPEN, UNDAMAGED PACKAGES

**CONTENTS:** MODEL NO. TFS3 - N- 3 WP  
 10 - Single use, disposable 15' Irrigation / aspiration tube set with probe  
 INDICATION FOR IRRIGATION  
 DURING LAPAROSCOPIC PROCEDURES

**CAUTION:** Federal (U.S.A.) Law restricts this device to sale by,  
 or on the order of a physician.  
 Recommended for Laparoscopic Irrigation  
 In adult and pediatric surgery.

THIS PRODUCT TO ONLY BE USED FOR  
 LAPAROSCOPIC IRRIGATION. DO NOT  
 USE FOR HYSTEROSCOPIC DISTENTION.

Distributed by:



**SANESE MEDICAL**  
 CORPORATION

888 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-8946 FAX (614) 291-4288

LABEL SIZE - 3.25 X 5.50  
 UL PDQ2 SPECIFICATIONS  
 TEXT COLOR- BLACK  
 BACKGROUND COLOR- WHITE

ECN	REV	DESCRIPTION	DATE	DWN
SANESE MEDICAL		TITLE	TFS3 - N- 3 WP CASE LABEL	
		PART NUMBER	10024	DWN DATE
SHEET 1 OF 1	SCALE	FULL	FILE	TFS3NWP_CDR MDI 12/12/94

101

	<p><b>READ THE OPERATIONS MANUAL BEFORE OPERATING THIS UNIT. FOR SALE BY, OR ON THE ORDER OF A HEALTH CARE PROFESSIONAL</b></p>
	<p><b>THIS UNIT CAN NOT BE OPERATED IN THE PRESENCE OF AN EXPLOSIVE ANAESTHETIC MIXTURE WITH AIR, OXYGEN OR NITROUS OXIDE.</b></p>

LABEL SIZE 1.90 X 3.00  
 CORNER RADI - .100  
 UL PDQ2 SPECIFICATIONS  
 TEXT COLOR- BLACK  
 BACKGROUND COLOR- WHITE  
 RED COLOR- PANTONE 032 RED  
 YELLOW COLOR - PANTONE YELLOW CV

ECN	REV	DESCRIPTION	DATE	DWN
<b>SANESE MEDICAL</b>		TITLE	AP WARNING LABEL	
		PART NUMBER	10012	DWN
		FILE	APWARN.CRD	MDT
SHEET 1 OF 1	SCALE	FULL		DATE
				12/12/94

100



This device is equipped with a "Hospital Grade" plug. Grounding reliability can only be achieved when the equipment is connected to an equivalent receptacle marked "Hospital Only" or "Hospital Grade".




### POWER REQUIREMENTS

INPUT - 110 VAC 60 HZ 1.0 AMP  
ONLY QUALIFIED PERSONNEL  
SHOULD PERFORM REPAIR  
TO THIS UNIT.

LABEL SIZE 2.00 x 3.60  
CORNER RADII - .100  
UL PDQ2 SPECIFICATIONS  
TEXT COLOR- BLACK  
BACKGROUND COLOR- WHITE  
YELLOW COLOR - PANTONE YELLOW CV

ECN	REV	DESCRIPTION	DATE	DWN		
SANESE MEDICAL		TITLE	POWER CORD WARNING LABEL			
		PART NUMBER	10013	DWN	DATE	
SHEET 1 OF 1	SCALE	FULL	FILE	POWER.CRD	MDI	12/12/94

103

	UL listed 17J5
SERIAL NO.	<input type="text"/>
MODEL NO.	<input type="text"/>
	
AIR	HEATER
For proper operation plug the pressure enclosure air line and the heater cable into the connectors.	

LABEL SIZE - 2.40 X 2.60  
 CORNER RADII - .100  
 UL PDQ2 SPECIFICATIONS  
 TEXT COLOR- BLACK  
 BACKGROUND COLOR- WHITE  
 RED - PANTONE S 73 -1

ECN	REV	DESCRIPTION	DATE	DWN
SANESE MEDICAL		TITLE CONNECTIONS LABEL		
		PART NUMBER 10025	DWN	DATE
SHEET 1 OF 1	SCALE FULL	FILE CONN.CDR	MDT	12/12/94

JDY

**PRESSURE**  
mm Hg

**POWER**

**THERMO-FLO**  
SYSTEMS

**HEATER**

**ON**

**LOW**

**LEVEL**

**RUN**

**LOAD**

**LABEL SIZE 2**  
OFF-SET COLOR PROCESS  
PDQ2 NOT REQUIRED  
TEXT COLOR- BLACK  
GREEN SHADES TBD PANTONE COLORS

ECH	REV	DESCRIPTION	DATE	DWN
<b>TITLE FRONT PANEL LABEL</b>				
PART NUMBER			10014	DWN
FILE			FRONT.CRD	MDT
SHEET 1 OF 1		SCALE	FULL	DATE
				12/12/94



# Safety Precautions

## Symbol Explanations



**READ THE OPERATIONS MANUAL BEFORE OPERATING THIS UNIT.  
FOR SALE BY, OR ON THE ORDER OF A HEALTH CARE PROFESSIONAL.**



**THIS UNIT CAN NOT BE OPERATED IN THE PRESENCE  
OF AN EXPLOSIVE ANAESTHETIC MIXTURE WITH AIR, OXYGEN  
OR NITROUS OXIDE.**



**Only qualified personnel should perform repair to this unit.  
This unit requires a power input of - 110 VAC 60 HZ 1.0 AMP  
This device is equipped with a "Hospital Grade" plug. Grounding  
reliability can only be achieved when the equipment is connected  
to an equivalent receptacle marked "Hospital Only" or "Hospital Grade".**

**Use only with TFS 3 Irrigation tube sets. These sets are equipped with  
lock and pierce components for the safe operation of the device. Never  
use parts from another manufacturer's unit even though they appear to fit.**

**Never service device while it is plugged into an outlet.**

**Never service device near water.**

**This device is designed to be operated from a power  
indicated on the marking label.**

**Do not overload wall outlets or use extension cords.**


**Do not allow anything to rest on the power cord or locate  
this product where the cord will be abused by persons or  
equipment moving on or near the cord.**

**This device is a Class II, Type B device.**

126



The Thermo-Flo System 3 was developed to eliminate the major problems faced everyday by the surgeons and O.R. unit Managers. These are cross contamination, inadvertent patient hypothermia, inadequate irrigation solution capacity, ease of set up, tube sets and disposable probes that leak, noise and physical space. Sanese Medical answers these problems with the Thermo-Flo System 3.



1) **CROSS CONTAMINATION** occurs when non sterile fluid enters the patient's abdominal cavity through the irrigation line. The Thermo - Flo System 3 eliminates cross contamination.

2) **INADVERTENT PATIENT HYPOTHERMIA** is the result of the combination of gas and fluids ranging in temperature from 65 to 72 degrees farenhite being introduced into the patient. The Thermo - Flo System can maintian the solution temperature at 100 degress farenhite. Though this system alone cannot prevent hypothermia, it will dramatically effect the reduction of the core temperature due to existing conditions of surgery.

**3) SOLUTION CAPACITY / EASE OF SET UP**

The Thermol - Flo System 3 meets the liter needs of the todays and tomorrows laparoscopic procedures. This amount provides less disruption of the procedure, by not replacing the solution, cap and tube inserts that promotes the inadvertent contamination by the person replacing the solution as often.

**KEY ADVANTAGES OF THE THERMOL-FLO SYSTEM 3**

- \* Delivers a strong stream of irrigation fluid
- \* No capital equipment necessary
- \* Choice of disposable or reusable probe
- \* Simplified O.R. set up / LV.
- \* Avoids Equipment purchase
- \* Quiet operation
- \* 24 hour, 7 days a week service available

**THERMO-FLOW SYSTEM 3 PRODUCT LINE**

- 1) The Thermo-Flo System 3 Irigation Pump.....TFS3 - one per case
- 2) TFS3-3WP Tube Set with Dispoable Probe.....TFS3-3WP one per case
- 3) TFS3-3WOP Tube Set with out Probe .....TFS3-3WOP one per case
- 4) TFS3-N-3WP Tube Set with Disposable Probe.....TFS3-N-3WP 10 per case
- 5) TFS3-n-3WP Tube set with out Probe.....TFS3-N-3WOP 10 per case

All Tube Sets are 15' with Aspiration line included

**SANESE MEDICAL CORPORATION**

222 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 271-2222 FAX (614) 271-4222

## **SPECIFICATIONS Model TFS3 \***

### **Power**

120 v ac 1 AMP or 220 V AC .5 AMP 50/60 Hz  
(95 to 250 V AC to 63 Hz)

### **Weight**

7 lbs. (with out stand)

### **Pressure Range**

100 mm Hg to 800 mm Hg + or - 10 mm Hg

### **Heater Temperature**

The heater is preset to 100 F° + or - 2 F°

Safety Shutdown Temperature- 105 F°.

### **Maximum Leakage Current**

Control Unit - 100 micro amps. (per UL 544)

Pressure Enclosure - 10 micro amps. (per UL 544)

Maximum Ground Impedance - 0.1 ohm ( per UL 544)

Fuse : Slow Blow 2.0 amp. ( Located on the back of control unit)

### **Environmental Operating Conditions**

40 F° to 150 F°, Humidity non-condensing

\* U. S. and International patent pending  
Specification subject to change without notice.

128

## DESCRIPTION OF OPERATION

The controller is plugged into the wall outlet. The power switch is turned on and the LOAD/RUN switch is placed in the LOAD position. A fresh preheated bag of irrigation fluid is placed in the pressure enclosure.

The irrigation bag is connected to the tubing set/probe through the stick fittings which protrude through a opening in the bottom of the pressure enclosure. The LOAD/RUN switch is placed in the run position.

When the power is first turned on the pressure setpoint will need to be set. The heater automatically comes on and goes to a set temperature of 100 F°. The heater is design to maintain the temperature of the irrigation fluid.

The user then depresses the fluid control valve on the probe to obtain irrigation fluid. After a time the FLUID LOW light will come on indicating that the level is low in the irrigation fluid bag.

When the bag becomes empty move the LOAD/RUN switch load position and replace the empty bag with a fresh bag of irrigation fluid.


The machine will remain on the last settings until they are changed by the user.

129

IRRIGATION PUMP SYSTEM

**THERMO-FLO**

SYSTEM 3



**SANESE MEDICAL**  
CORPORATION

The new Thermo-Flo System 3 was developed to eliminate the major problems faced everyday by the surgeon and the O.R. unit manger.

150

# **OPERATING THE TFS3 IRRIGATION SYSTEM**

## **DESCRIPTION OF THE CONTROLS**

**Power Switch:** The power switch is located on the front panel. It controls the power to unit.

**Pressure Adjust Switch:** Located to the right of the displays a three position rocker switch which is used to adjust the pressure. The center position is off. The top position (up arrow) and bottom position (down arrow) are momentary. Pressing the UP arrow increases the pressure. Pressing the bottom (down arrow) will decrease the pressure.

**Load/Run Switch:** This switch is a two position rocker. When the switch is set to the Load position the device will drop the air pressure in the Pressure Enclosure so the bag can be reloaded. When the switch is in the Run position the Pressure Enclosure will inflate to the set point on the display.

**Display :** The display is a 16 character single line LCD.

**Low Fluid LED:** This is a panel light on the front panel that indicates when the fluid level is low and prompts the user to prepare a new bag of solution.

**Heater On LED:** This is a panel light that indicates when the heater is on. If this light is not on the heater is not maintaining the Pre-set temperature of 100° F.

**Air Connector:** On the Back of panel of the unit a Colder Connector which plugs into a mating connector on the Pressure Enclosure.

**Heater Electrical Connector:** On the back panel of the controller is located a connector used to electrically connect the heater and the temperature sensor to the controller.

151

Food and Drug Administration  
Center for Devices and  
Radiological Health  
Office of Device Evaluation  
Document Mail Center (HFZ-401)  
9200 Corporate Blvd.  
Rockville, Maryland 20850

December 15, 1994

SANESE MEDICAL CORPORATION  
885 NORTHWEST BLVD.  
COLUMBUS, OH 43212  
ATTN: CHRISTOPHER N. SANESE

510(k) Number: K941507  
Product: THEMO-FLO SYSTEM  
3  
LAPAROSCOPIC/HYS  
TEROSCOPIC

Extended Until: 17-FEB-95

Based on your recent request, an extension of time has been granted for you to submit the additional information we requested.

If the additional information is not received by the "Extended Until" date shown above your premarket notification will be considered withdrawn.

If you have procedural or policy questions, please contact the Division of Small Manufacturers Assistance at (301) 443-6597 or at their toll-free number (800) 638-2041, or contact me at (301) 594-1190.

Sincerely yours,

Marjorie Shulman  
Supervisory Consumer Safety Officer  
Premarket Notification Section  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

132



# SANESE MEDICAL CORPORATION

MANUFACTURER OF MEDICAL PRODUCTS

RECEIVED

13 DEC 94 14 03

FDA/CDRH/OBE/DHC

12-9-94

Colin M. Pollard  
Food and Drug Administration  
Center for Devices and  
Radiological Health  
Document Mail Center (HFZ-401)  
9200 Corporate Boulevard  
Rockville, Maryland 20850

RE:K941507 / A  
Thermo - Flo System 3  
Dated: July 5, 1994  
Received: July 6

Dear Mr. Pollard,

This letter is to formally request a 60 day extension of our December 17, delivery compliance as requested in your letter of November 17, 1994.

This letter is a follow up to the verbal request made of Donna -Bea Tillman in our phone conversation of December 5, 1994.

Thank you very much for your consideration of this request.

Sincerely,

Christopher N. Sanese  
Regulatory Correspondent  
SANESE MEDICAL CORPORATION

ENC: 4 pages

CNS/III

CC: John Jevicky, Esq.  
DINSMORE & SHOHL



133





Food and Drug Administration  
1390 Piccard Drive  
Rockville MD 20850

Mr. Christopher N. Sanese  
Regulatory Correspondent  
Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, Ohio 43212

NOV 17 1994

Re: K941507/A  
Thermo-Flo System 3  
Dated: July 5, 1994  
Received: July 6, 1994

Dear Mr. Sanese:

We have reviewed your Section 510(k) notification of intent to market the device referenced above. We cannot determine if the device is substantially equivalent to a device marketed prior to May 28, 1976, the enactment date of the Medical Device Amendments, based solely on the information you provided.

[REDACTED] (b)(4) [REDACTED] [REDACTED] [REDACTED]

(b)(4) Trade Secret Process - Product Specification

7521

(b)(4) Trade Secret Process - Product Specification



135

(b)(4) Trade Secret Process - Product Specification



136

Page 4 - Mr. Christopher N. Sanese

bcc: HFZ-401 DMC  
HFZ-404 510(k) Staff  
HFZ- Division  
D.O.

DBTillman  
C:\WP51\FILES\K941507.AI  
Draft 11-15-94 ens

137

Mr. Christopher N. Sanese  
Regulatory Correspondent  
Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, Ohio 43212

Re: K941507/A  
Thermo-Flo System 3  
Dated: July 5, 1994  
Received: July 6, 1994

Dear Mr. Sanese:

We have reviewed your Section 510(k) notification of intent to market the device referenced above. We cannot determine if the device is substantially equivalent to a device marketed prior to May 28, 1976, the enactment date of the Medical Device Amendments, based solely on the information you provided.

(b)(4) Trade Secret Process - Product Specification [REDACTED]:

(b)(4) Trade Secret Process - Product Specification [REDACTED]

138

(b)(4) Trade Secret Process - Product Specification



1309

(b)(4) Trade Secret Process - Product Specification



140

DEPARTMENT OF HEALTH AND HUMAN SERVICES

Page 4 - Mr. Christopher N. Sanese

bcc: HFZ-401 DMC  
HFZ-404 510(k) Staff  
HFZ- Division  
D.O.

DBTillman  
C:\WP51\FILES\K941507.AI  
Draft 11-15-94 ens

F" E  
COPY

141

OFFICE	SURNAME	DATE	OFFICE	SURNAME	DATE	OFFICE	SURNAME	DATE
Z-470	Tillman	11/5/94						
Z470	Palard	11/16/94						





# Memorandum

Date: November 15, 1994

From: REVIEWER(S) - NAME(S) Donna-Bea Tillman

Subject: 510(k) NOTIFICATION K94/507/S1

To: THE RECORD

It is my recommendation that the subject 510(k) Notification:

- (A) Is substantially equivalent to marketed devices.
- (B) Requires premarket approval. NOT substantially equivalent to marketed devices.
- (C) Requires more data.
- (D) Other (e.g., exempt by regulation, not a device, duplicate, etc.)

Additional Comments:

Is this device subject to Postmarket Surveillance? Yes  No

This 510(k) contains: (check appropriate box(es))

- A 510(k) summary of safety and effectiveness, or
- A 510(k) statement that safety and effectiveness information will be made available
- The required certification and summary for class III devices

The submitter requests under 21 CFR 807.95:\*

- No Confidentiality
- Confidentiality for 90 days
- Continued Confidentiality exceeding 90 days

Predicate Product Code w/panel and class:

Additional Product Code(s) w/Panel (optional):

REVIEW: Colin M Pollard  
(BRANCH CHIEF)


06D13 11/16/94  
BRANCH CODE (DATE)

FINAL REVIEW: \_\_\_\_\_  
(DIVISION DIRECTOR) (DATE)

142

\*DOES NOT APPLY TO ANY "SE" DECISIONS

# 510(k) "SUBSTANTIAL EQUIVALENCE" DECISION-MAKING PROCESS (DETAILED)



\* 510(k) submissions compare new devices to marketed devices. FDA requests additional information if the relationship between marketed and "predicate" (pre-Amendments or reclassified post-Amendments) devices is unclear.

\*\* This decision is normally based on descriptive information alone, but limited testing information is sometimes required.

\*\*\* Data may be in the 510(k), other 510(k)s, the Center's classification files, or the literature.

# OGDB

## K941507/A

**Reviewer:** Donna-Bea Tillman  
Biomedical Engineer

**Division/Branch:** DRAERD/ADOU/OGDB  
(HFZ-470)

**Proprietary Trade Name:** Thermo-Flo System 3  
**Common Name:** Laparoscopic Irrigator

**Product to which compared:** Cabot Niagra High-Flow Irrigator (K924530)

**Applicant:** Sanese Medical Corporation  
885 Northwest Blvd  
Columbus, OH 43212

**Contact:** Christopher N. Sanese, Regulatory Correspondent  
**Phone:** (614) 291-8848

## DEVICE DESCRIPTION

1. *Intended Use:*

The Thermo-Flo System 3 is intended for irrigation/aspiration during laparoscopic surgery (Note: In a October 15, 1994, amendment, the firm withdrew its claims for hysteroscopy).

2. *Physical Description:*

	YES	NO
· Is the device life-supporting or life sustaining?	—	<u>X</u>
· Is the device implanted (short-term or long-term)?	—	<u>X</u>
· Does the device design use software?	<u>X</u>	—
· Is the device sterile?	—	<u>X</u>
· Is the device single use?	—	<u>X</u>
· Is the device home use?	—	<u>X</u>
· Is the device for prescription?	<u>X</u>	—
· Does the device contain a drug or biological product as a component?	—	<u>X</u>
· Is this device a kit?	—	<u>X</u>

144

(b)(4) Trade Secret Process - Product Specification



## REVIEW ANALYSIS

(b)(4) Trade Secret Process - Product Specification



145

(b)(4) Trade Secret Process - Product Specification



1-10

(b)(4) Trade Secret Process - Product Specification



117

· Labeling

The Operator's manual does not contain indications for use, contraindications, or warnings. In addition, the cleaning instructions: "Do not use abrasive cleaning and disinfecting solutions" are inadequate.

- x. Please modify your labeling to include specific indications for use, contraindications, warnings, and detailed cleaning instructions.

## REVIEWER RECOMMENDATION

Hold for Additional Information

Donna-Bea Tillman      11/15/94  
Donna-Bea Tillman      Date

Colin M. Pollard      11/16/94  
Colin M. Pollard      Date  
Chief, Ob/Gyn Devices Branch

Concur  
 Do not concur.  
Comments:

1148

DIVISION OF REPRODUCTIVE, ABDOMINAL, EAR, NOSE AND THROAT  
AND RADIOLOGICAL DEVICES

MEMORANDUM OF TELEPHONE CONVERSATION

Date: October 11, 1994

BETWEEN: Donna-Bea Tillman

AND: Christopher Sanese

TITLE: Regulatory Agent

COMPANY: Sanese Medical

PHONE NUMBER: 614-291-8848

FAX NUMBER: \_\_\_\_\_

DOCUMENT NUMBER: K941507

SUMMARY:

(b)(4) Trade Secret Process - Product Specification



Signed: Donna-Bea Tillman

Original to: K941507

Copy to: dbtlog

149



K941507/A2

**SANESE MEDICAL**  
CORPORATION

MANUFACTURER OF MEDICAL PRODUCTS

RECEIVED

19 OCT 94 14 33

FDA/CDRH/ODE/DMC

*Mr. Collin Pallard/ Ms. Donna B. Tillman  
Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center (HFZ-401)  
1390 Piccard Drive  
Rockville, Maryland 20850*

October 15, 1994

**RE: K941507\_**  
*Thermo-Flo System 3  
Laparoscopic/Hysteroscopic  
Irrigation System and  
Disposable Irrigation/  
Aspiration Probe and Tubing  
Dated: March 18, 1994*

VIA FAX & FEDERAL EXPRESS

Dear Mr. Pallard & Ms. Tillman

(b)(4) Trade Secret Process - Product Specification



Respectfully,

Christopher N. Sanese  
REGULATORY AGENT

ISO

LOT NO. 0300



**CONTENTS:**

1 - MODEL TFS3 LAPAROSCOPIC IRRIGATION SYSTEM

INDICATED FOR IRRIGATION DURING LAPAROSCOPIC PROCEDURES

**Distributed by:**

**SANESE MEDICAL**  
CORPORATION

885 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-8848 FAX (614) 291-1389

**CAUTION:** Federal ( U.S.A. ) Law restricts this device to the sale by,  
or on the order of a physician.

**NON - STERILE: RECOMMENDED FOR LAPAROSCOPIC IRRIGATION  
IN ADULT AND PEDIATRIC SURGERY**

151

OPERATIONS MANUAL



LAPAROSCOPIC IRRIGATION SYSTEM

152

**SANESE MEDICAL**  
CORPORATION

885 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-8848 FAX (614) 291-1389

## ***SAFETY PRECAUTIONS***

**\*\* USE ONLY WITH TFS 3 IRRIGATION TUBE SETS \*\***

- \* Do not use in volatile atmosphere.
- \* Use only with 120 volt grounded receptacle.
- \* Do not immerse unit.

## **READ BEFORE ADJUSTING AND OPERATING UNIT**

It is recommended that all O.R. Unit Managers and others responsible for the operation of this unit become thoroughly familiar with its capabilities and proper operating procedures prior to actual patient use.

## ***PRODUCT DESCRIPTION***

The **THERMO-FLO SYSTEM 3™** LAPAROSCOPIC IRRIGATION SYSTEM consists of the following components:

- 1 - **THERMO - FLO SYSTEM 3™** control device
- 1 - **THERMO - FLO SYSTEM 3™** pressure enclosure device
- 2 - Transportation pouch

\*U. S. and International patent pending  
**SANESE MEDICAL**  
CORPORATION

ISS

CONTENTS: **MODEL NO. TFS3 - N 2WP**

1 - Single use, disposable 15' irrigation / aspiration tube set with probe

**LOT NO. 0500**

**FOR SINGLE USE ONLY**

**CONTENTS STERILE IN UNOPENED, UNDAMAGED PACKAGES.**

**INDICATED FOR IRRIGATION DURING LAPAROSCOPIC PROCEDURES**

***Distributed by:*** **SANESE MEDICAL**  
CORPORATION

885 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-8848 FAX (614) 291-1389

**RECOMENDED FOR ADULT AND PEDIATRIC LAPAROSCOPIC SURGERY**

***CAUTION:*** *Federal ( U.S.A. ) law restricts this device to sale by, or on the order of a physician.*

**THIS PRODUCT TO ONLY BE USED FOR LAPAROSCOPIC IRRIGATION.  
DO NOT USE FOR HYSTEROSCOPIC DISTENTION.**

154

**CONTENTS: MODEL NO. TFS3 - 3 WOP**

1 - Single use, disposable 15' irrigation / aspiration tube set.

**LOT NO. 0500**

**FOR SINGLE USE ONLY**

**CONTENTS STERILE IN UNOPENED, UNDAMAGED PACKAGES.**

**INDICATED FOR IRRIGATION DURING LAPAROSCOPIC PROCEDURES**

*Distributed by:* **SANESE MEDICAL**  
CORPORATION

885 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-8848 FAX (614) 291-1389

**RECOMENDED FOR ADULT AND PEDIATRIC LAPAROSCOPIC SURGERY**

***CAUTION:** Federal ( U.S.A.) law restricts this device to sale by, or on the order of a physician.*

**THIS PRODUCT TO ONLY BE USED FOR LAPAROSCOPIC IRRIGATION.  
DO NOT USE FOR HYSTEROSCOPIC DISTENTION.**

155

**CONTENTS: MODEL NO. TFS3 - 3 WP**

1 - Single use, disposable 15' irrigation / aspiration tube set with probe

**LOT NO. 0500**

**FOR SINGLE USE ONLY**

**CONTENTS STERILE IN UNOPENED, UNDAMAGED PACKAGES.**

**INDICATED FOR IRRIGATION DURING LAPAROSCOPIC PROCEDURES**

*Distributed by:*

**SANESE MEDICAL**  
CORPORATION

885 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-8848 FAX (614) 291-1389

**RECOMENDED FOR ADULT AND PEDIATRIC LAPAROSCOPIC SURGERY**

*CAUTION: Federal ( U.S.A.) law restricts this device to sale by, or on the order of a physician.*

**THIS PRODUCT TO ONLY BE USED FOR LAPAROSCOPIC IRRIGATION.  
DO NOT USE FOR HYSTEROSCOPIC DISTENTION.**

156

USAGE RECOMMENDATION

The **THERMO-FLO SYSTEM 3™** was designed for Laparoscopic surgery and is only to be used by licensed, certified surgical personnel.

RECOMMENDED CARE

Do not leave unit plugged in.

Do not use abrasive cleaning and disinfecting solutions.

LIMITED WARRANTY

Sanese Medical Corporation ("Sanese") warrants the **THERMO-FLO SYSTEM 3™** to be free from defects in material and workmanship for a period of one (1) year from date of initial purchase. No warranty shall apply if the goods have been damaged by accident, abuse, misuse, or misapplication, or as a result of service or modification by other than a person authorized by Sanese. Sanese's liability under this warranty, and buyer's exclusive remedy is limited to the cost of materials and labor to repair defective goods, or to their replacement, and in no event shall exceed the purchase price. Repair or replacement of defective goods under warranty will be made only upon their return to Sanese after notice to Sanese and Buyer's receipt of shipping instructions. Upon receipt of goods returned under warranty, Sanese will inspect the goods, will notify buyer of the extent of repair or replacement which Sanese will perform under warranty, which shall be conclusive of Sanese's liability.

Sanese is not responsible for incidental consequential damages resulting from the breach of any expressed or implied warranty, including damage to property and to the extent permitted by law, damages for personal injury. The warranty contained herein is in lieu of all warranties, expressed or implied. No statement of any representative shall extend Sanese's liability as herein established or limited.

157

**SANESE MEDICAL**  
CORPORATION



# SANESE MEDICAL CORPORATION

MANUFACTURER OF MEDICAL PRODUCTS

## FACSIMILE TRANSMISSION

301-594-2359

9 pages

### ATTENTION:

*Mr. Collin Pallard / Ms. Donna B. Tillman  
Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center (HFZ-401)  
1390 Piccard Drive  
Rockville, Maryland 20850*

*October 15, 1994*

*510(k) Number: K941507  
Received: 29-Mar-94  
Product:  
Thermo-Flo System 3  
Laparoscopic /Hysteroscopic*

Christopher N. Sancec  
Regulatory Correspondent

158

# SANESE MEDICAL CORPORATION

MANUFACTURER OF MEDICAL PRODUCTS

*Mr. Collin Pallard/Ms. Donna B. Tillman  
Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center (HFZ-401)  
1390 Piccard Drive  
Rockville, Maryland 20850*

*October 15, 1994*

*RE: K941507  
Thermo-Flo System 3  
Laparoscopic/Hysteroscopic  
Irrigation System and  
Disposable Irrigation/  
Aspiration Probe and Tubing  
Dated: March 18, 1994*


**VIA FAX & FEDERAL EXPRESS**

Dear Mr. Pallard & Ms. Tillman

(b)(4) Trade Secret Process - Product Specification



Respectfully,



Christopher N. Sanese  
REGULATORY AGENT

159









**OPERATIONS MANUAL**



**LAPAROSCOPIC IRRIGATION SYSTEM**

161

## ***SAFETY PRECAUTIONS***

**\*\* USE ONLY WITH TFS 3 IRRIGATION TUBE SETS \*\***

- \* Do not use in volatile atmosphere.
- \* Use only with 120 volt grounded receptacle.
- \* Do not immerse unit.

## **READ BEFORE ADJUSTING AND OPERATING UNIT**

It is recommended that all O.R. Unit Managers and others responsible for the operation of this unit become thoroughly familiar with its capabilities and proper operating procedures prior to actual patient use.

## ***PRODUCT DESCRIPTION***

The **THERMO-FLO SYSTEM 3™** LAPAROSCOPIC IRRIGATION SYSTEM consists of the following components:

- 1 - **THERMO - FLO SYSTEM 3™** control device
- 1 - **THERMO - FLO SYSTEM 3™** pressure enclosure device
- 2 - Transportation pouch

U. S. and International patent pending.  
**SANESE MEDICAL**  
CORPORATION

685 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-8848 FAX (614) 291-1389

165



### USAGE RECOMMENDATION

The **THERMO-FLO SYSTEM** 3™ was designed for Laparoscopic surgery and is only to be used by licensed, certified surgical personnel.

### RECOMMENDED CARE

Do not leave unit plugged in.

Do not use abrasive cleaning and disinfecting solutions.

### LIMITED WARRANTY

Sanese Medical Corporation ("Sanese") warrants the **THERMO-FLO SYSTEM** 3™ to be free from defects in material and workmanship for a period of one (1) year from date of initial purchase. No warranty shall apply if the goods have been damaged by accident, abuse, misuse, or misapplication, or as a result of service or modification by other than a person authorized by Sanese. Sanese's liability under this warranty, and buyer's exclusive remedy is limited to the cost of materials and labor to repair defective goods, or to their replacement, and in no event shall exceed the purchase price. Repair or replacement of defective goods under warranty will be made only upon their return to Sanese after notice to Sanese and Buyer's receipt of shipping instructions. Upon receipt of goods returned under warranty, Sanese will inspect the goods, will notify buyer of the extent of repair or replacement which Sanese will perform under warranty, which shall be exclusive of Sanese's liability.

Sanese is not responsible for incidental consequential damages resulting from the breach of any expressed or implied warranty, including damage to property and to the extent permitted by law, damages for personal injury. The warranty contained herein is in lieu of all warranties, expressed or implied. No statement of any representative shall extend Sanese's liability as herein established or limited.

**SANESE MEDICAL**  
CORPORATION

Hardcopy

**SANESE MEDICAL**  
CORPORATION

MANUFACTURER OF MEDICAL PRODUCTS

K941507/A1

Mr. Collin Pallard  
Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center (HFZ-401)  
1390 Piccard Drive  
Rockville, Maryland 20850

August 10, 1994

RE: K941507  
Thermo-Flo System 3  
Laparoscopic/Hysteroscopic  
Irrigation System and  
Disposable Irrigation/  
Aspiration Probe and Tubing  
Dated: March 18, 1994

FDA/CDRH/CDE/DNC

15 AUG 94 09 03

RECEIVED

VIA FAX & REGULAR MAIL


Dear Mr. Pallard,

Here is the additional information you requested regarding the distal tip of the probe.

The material for the distal tip of the probe will be AISI stainless steel type 304 (for disposable) & 316 (for reusable). The Tube set with and with out probe will be sterilized via gamma radiation as determined by AAMI method 1 or 2.5 MRads, which ever is greater, with a Sterility Assurance level of 10<sup>-6</sup>.

Thank you for your time and efforts, it is greatly appreciated.

Sincerely,



Christopher N. Sanese  
Regulatory Correspondent

162

Food and Drug Administration  
Center for Devices and  
Radiological Health  
Office of Device Evaluation  
Document Mail Center (HFZ-401)  
1390 Piccard Drive  
Rockville, Maryland 20850

July 07, 1994

SANESE MEDICAL CORPORATION  
885 NORTHWEST BLVD.  
COLUMBUS, OH 43212  
ATTN: CHRISTOPHER N. SANESE

510(k) Number: K941507  
Product: THEMO-FLO SYSTEM  
3  
LAPAROSCOPIC/HYS  
TEROSCOPIC

The additional information you have submitted has been received.

We will notify you when the processing of this submission has been completed or if any additional information is required. Please remember that all correspondence concerning your submission MUST be sent to the Document Mail Center (HFZ-401) at the above letterhead address. Correspondence sent to any address other than the one above will not be considered as part of your official premarket notification submission. Because of equipment and personnel limitations we cannot accept telefaxed material as part of your official premarket notification submission, unless specifically requested of you by an FDA official.

The Safe Medical Devices Act of 1990, signed on November 28, states that you may not place this device into commercial distribution until you receive a letter from FDA allowing you to do so. As in the past, we intend to complete our review as quickly as possible. Generally we do so 90 days. However, the complexity of a submission or a requirement for additional information may occasionally cause the review to extend beyond 90 days. Thus, if you have not received a written decision or been contacted within 90 days of our receipt date you may want to check with FDA to determine the status of your submission.

If you have procedural or policy questions, please contact the Division of Small Manufacturers Assistance at (301) 443-6597 or at their toll-free number (800) 638-2041, or contact me at (301) 594-1190.

Sincerely yours,

Marjorie Shulman  
Supervisory Consumer Safety Officer  
Premarket Notification Section  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

108

**SANESE MEDICAL**  
CORPORATION

MANUFACTURER OF MEDICAL PRODUCTS

*Mr. Collin Pallard  
Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center (HFZ-401)  
1390 Piccard Drive  
Rockville, Maryland 20850*

July 5, 1994


*RE: K941507  
Thermo-Flo System 3  
Laparoscopic/Hysteroscopic  
Irrigation System and  
Disposable Irrigation/  
Aspiration Probe and Tubing  
Dated: March 18, 1994*

Dear Mr. Pallard,

(b)(4) Trade Secret Process - Product Specification



Sincerely,

  
Christopher N. Sanese  
Regulatory Correspondent

*K941507/S*

**SANESE MEDICAL CONTROLLED DOCUMENT**  
**Title: Verification/Validation description - System 3 pump**  
**Written: K.A. Kelly**  
**Last Revised on ECN: N/A**

**No: eo001-03Rev A**  
**Page: 1 of 2**  
**Releasd: 6/28/94**

b(4)Trade Secret Process-Product Specs



170

**SANESE MEDICAL CONTROLLED DOCUMENT**  
**Title: Verification/Validation description - System 3 pump**  
**Written: K.A. Kelly**  
**Last Revised on ECN: N/A**

**No: eo001-03Rev A**  
**Page: 2 of 2**  
**Releasd: 6/28/94**

b(4)Trade Secret Process-Product Specs



171

**SANESE MEDICAL CONTROLLED DOCUMENT**

**Title: Hazard Analysis - System 3 Irrigation Pump**

**Written: K.A. Kelly**

**Last Revised on ECN: N/A**

**No: ci001-02Rev A**

**Page: 1 of 5**

**Release Date: 6/20/94**

b(4)Trade Secret Process-Product Specs



172

**SANESE MEDICAL CONTROLLED DOCUMENT**

**Title: Hazard Analysis - System 3 Irrigation Pump**

**Written: K.A. Kelly**

**Last Revised on ECN: N/A**

**No: ci001-02Rev A**

**Page: 2 of 5**

**Release Date: 6/20/94**

b(4)Trade Secret Process-Product Specs



173



**SANESE MEDICAL CONTROLLED DOCUMENT**

**Title: Hazard Analysis - System 3 Irrigation Pump**

**Written: K.A. Kelly**

**Last Revised on ECN: N/A**

**No: ei001-02Rev A**

**Page: 3 of 5**

**Release Date: 6/20/94**

b(4)Trade Secret Process-Product Specs



174

**SANESE MEDICAL CONTROLLED DOCUMENT**

**Title: Hazard Analysis - System 3 Irrigation Pump**

**Written: K.A. Kelly**

**Last Revised on ECN: N/A**

**No: ei001-02Rev A**

**Page: 4 of 5**

**Release Date: 6/20/94**

b(4)Trade Secret Process-Product Specs



175

**SANESE MEDICAL CONTROLLED DOCUMENT**  
**Title: Hazard Analysis - System 3 Irrigation Pump**  
**Written: K.A. Kelly**  
**Last Revised on ECN: N/A**

**No: ei001-02Rev A**  
**Page: 5 of 5**  
**Release Date: 6/20/94**

b(4)Trade Secret Process-Product Specs



176



Food and Drug Administration  
1390 Piccard Drive  
Rockville, MD 20850

JUN - 7 1994

Mr. Christopher N. Sanese  
Regulatory Correspondent  
Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, Ohio 43212

Re: K941507  
Thermo-Flo System 3, Laparoscopic/  
Hysteroscopic Irrigation System and  
Disposable Irrigation/Aspiration  
Probe and Tubing  
Dated: March 18, 1994  
Received: March 29, 1994

Dear Mr. Sanese:

We have conducted an administrative review of your Section 510(k) notification of intent to market the device referenced above. This administrative review is part of a pilot screening program in the Division of Reproductive, Abdominal, Ear, Nose and Throat, and Radiological Devices (DRAERD). Our review indicates that your 510(k) is administratively incomplete and we are placing your 510(k) on hold. We believe that basic information is necessary for us to begin our technical review and to determine whether or not this device is substantially equivalent to a legally marketed device with regard to its safety and effectiveness. In order for us to begin the technical review of your submission, we require the information indicated on the enclosed DRAERD Premarket Notification 510(k) Reviewer's Screening Checklist. An explanation of each of the items on the Checklist is enclosed. Also enclosed for your information is a copy of the Office of Device Evaluation's (ODE) "510(k) 'Substantial Equivalence' Decision-Making Process (Detailed)" flow-chart, and a copy of the DRAERD Draft Guidance for the Content of Premarket Notifications. Since your 510(k) premarket notification has not been technically reviewed, additional information may be required during the review process and the file may again be placed on hold.

The additional information should be submitted in duplicate, referencing the 510(k) number above to:


Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center (HFZ-401)  
1390 Piccard Drive  
Rockville, Maryland 20850

You may not market this device until you have provided adequate information as required by 21 CFR 807.87(f) and (h), and you have received a letter from FDA allowing you to do so. If you market the device without conforming to these requirements, you will be in violation of the Federal Food, Drug, and Cosmetic Act (Act). You may, however, distribute this device for investigational purposes to obtain clinical data if needed to establish substantial equivalence. Clinical investigations of this device must be conducted in accordance with the investigational device exemptions (IDE) regulations.

If the requested information is not received within 30 days, we will consider your premarket notification to be withdrawn and your submission will be deleted from our system. If you submit the requested information after 30 days, it will be considered and processed as a new 510(k); therefore, all information previously submitted must be resubmitted so that your new 510(k) is complete.

If you have any questions concerning the contents of this letter, or comments/questions concerning this DRAERD pilot screening program, please contact Robert Gatling at (301) 594-1220. If you need information or assistance concerning the IDE or other regulations or copies of guidance documents, please contact the Division of Small Manufacturers Assistance at their toll free number 1-800-638-2041 or at (301) 443-6597.

Sincerely yours,



Lillian Yin, Ph.D.  
Director, Division of Reproductive,  
Abdominal, Ear, Nose and Throat,  
and Radiological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosures

178

510(K) ROUTE SLIP

510(k) NUMBER K941507 PANEL OB DIVISION DRAER BRANCH OGDB

TRADE NAME THEMO-FLO SYSTEM 3 LAPAROSCOPIC/HYSTEROSCOPIC SYSTEM

COMMON NAME LAPAROSCOPIC IRRIGATION PUMP SYSTEM

PRODUCT CODE \_\_\_\_\_

APPLICANT SANESE MEDICAL CORPORATION

SHORT NAME SANEMEDI

CONTACT CHRISTOPHER N SANESE

DIVISION \_\_\_\_\_

ADDRESS 885 NORTHWEST BLVD.

COLUMBUS, OH 43212

PHONE NO. (614) 291-8848

FAX NO. (\_\_\_\_) \_\_\_\_-\_\_\_\_

MANUFACTURER SANESE MEDICAL CORPORATION

REGISTRATION NO. \_\_\_\_\_

DATE ON SUBMISSION 16-MAR-94

DATE DUE TO 510(K) STAFF 12-JUN-94

DATE RECEIVED IN ODE 29-MAR-94

DATE DECISION DUE 27-JUN-94

DECISION \_\_\_\_\_

DECISION DATE \_\_\_\_\_

CORRESPONDENCE SENT DUE BACK

C001            06-JUN-94            06-JUL-94            HOLD LETTER

179

Memorandum

Date

June 2, 1994

From

REVIEWER(S) - NAME(S)

Donna-Bea Tillman

Subject

510(k) NOTIFICATION

K941507

To

THE RECORD

It is my recommendation that the subject 510(k) Notification:

- (A) Is substantially equivalent to marketed devices.
- (B) Requires premarket approval. NOT substantially equivalent to marketed devices.
- (C) Requires more data. Failed Checklist
- (D) Other (e.g., exempt by regulation, not a device, duplicate, etc.)

I have called the mfr. about this. comp

Additional Comments:

Is this device subject to Postmarket Surveillance? Yes  No

This 510(k) contains: (check appropriate box(es))

- A 510(k) summary of safety and effectiveness, or
- A 510(k) statement that safety and effectiveness information will be made available
- The required certification and summary for class III devices

The submitter requests under 21 CFR 807.95:\*

Predicate Product Code w/panel and class:

- No Confidentiality
- Confidentiality for 90 days
- Continued Confidentiality exceeding 90 days

Additional Product Code(s) w/Panel (optional):

REVIEW: Colin M. Pollard  
(BRANCH CHIEF)

09DB | 6/6/94  
BRANCH CODE (DATE)

FINAL REVIEW: Robert J. L. Linn  
(DIVISION DIRECTOR)

6/6/94  
(DATE)


\*DOES NOT APPLY TO ANY "SE" DECISIONS

Revised 11/18/91

180

# 510(k) "SUBSTANTIAL EQUIVALENCE" DECISION-MAKING PROCESS (DETAILED)

KAPLAN



- 510(k) submissions compare new devices to marketed devices. FDA requests additional information if the relationship between marketed and "predicate" (pre-Amendments or reclassified post-Amendments) devices is unclear.
- \*\* This decision is normally based on descriptive information alone, but limited testing information is sometimes required.
- \*\*\* Data may be in the 510(k), other 510(k)s, the Center's classification files, or the literature.



RRG/LLD 1/6/93  
Rev. 3/30/94

DRAERD Premarket Notification 510(k)  
Screening Checklist

510(k) Number &

Device Name K941507 Thermo-Flo System 3 Hysteroscopic/Laparoscopic Irrigation System and Disposable Irrigation/Aspiration Probe and Tubing  
Company Sanese Medical Corp.

ITEM	PRESENT		NEEDED (Y/N/?)
	Yes	No	
1. General information ( i.e., trade & classification name, Est. Reg. No., device class, meets special controls or a performance standards, etc.)	<u>X</u>	—	—
Reason for 510(k) - new device or modification	<u>X</u>	—	—
Identification of legally marketed equivalent device	<u>X</u>	—	—
2. Proposed Labeling, Labels, Advertisements	<u>X</u>	—	—
Description of new device/modification	<u>X</u>	—	—
Intended use statement	<u>X</u>	—	—
Diagrams, Engineering Drawings, Photographs	<u>X</u>	—	—
3. Comparison of similarities/differences to named legally marketed equivalent device	<u>X</u>	—	—
Equivalent Device Labeling, Labels, Advertising	<u>X</u>	—	—
Intended use of equivalent device	<u>X</u>	—	—
4. List of all patient contacting materials in new device	<u>X</u>	—	—
Comparison of materials to equivalent device	<u>X</u>	—	—
5. Biocompatibility information/data for <u>patient<sup>1</sup> contacting materials</u> , OR, Certification - identical material/formulation	—	<u>X</u>	<u>Y</u>
	—	<u>X</u>	—
6. Performance data: Bench data <sup>2</sup>	—	<u>X</u>	<u>Y</u>
Animal data	—	<u>X</u>	<u>N</u>
Clinical data	—	<u>X</u>	<u>?</u>
7. Sterilization information	<u>X</u>	—	—
8. Software validation & verification* *See attached guidance document	—	<u>X</u>	<u>Y</u>
9. 510(k) summary or statement	<u>X</u>	—	—
10. If Class III, Class III Certification & Summary	—	—	<u>NA</u>
11. If kit, kit certification	—	—	<u>NA</u>

182

<sup>1</sup> **Biocompatibility information**

(b)(4) Trade Secret Process - Product Specification

A large black rectangular redaction box covering the content of the first section.

<sup>2</sup> **Performance data**

(b)(4) Trade Secret Process - Product Specification

A large black rectangular redaction box covering the content of the second section.

183

FOR REVIEWER'S USE ONLY

RRG 9/24/93

DRAERD Premarket Notification 510(k)  
SUPPLEMENTAL  
Reviewer's Screening Checklist

DRAERD has been given the go ahead to continue with the DRAERD Premarket Notification 510(k) Screening Checklist program rather than switching to the ODE Premarket Notification (510(k)) Checklist for Acceptance Decision. However, some items appear in the ODE Checklist that were not in the early version of the DRAERD Checklist or Explanation of the Checklist. Therefore, the following items should be included as part of the DRAERD screening process:

510(k) Number: K941507 TIER (Circle) I / II / III

Expedited Review Requested: Y / N Granted: Y / N OR, FDA Identified Expedited: Y / N

ITEM	Yes	No
1. Is the product a device?	<u>X</u>	<u>  </u>
2. Is the device exempt from 510(k) by regulation or policy?	<u>  </u>	<u>X</u>
3. Are you aware that this device has been the subject of a previous NSE decision? If yes, does this new 510(k) address the NSE Issue(s) (e.g., performance data)?	<u>  </u>	<u>X</u>
4. Are you aware of the submitter being the subject of an integrity investigation? If yes, consult the ODE Integrity Officer, and has the ODE Integrity Officer given permission to proceed with the review?	<u>  </u>	<u>X</u>
5. Is there a specific guidance document for this device or device issue(s)?	<u>X</u>	<u>  </u>

In addition, the following item is new to the 510(k) review process. It will not be counted as a screening deficiency since it is new and "unknown" to the industry. It should be identified as a deficiency and requested as part of the technical review. The Explanation of the DRAERD Screening Checklist has been modified to include this information.

6. Address of manufacturing facility/facilities, and if applicable, sterilization site(s). X   

Administrative Reviewer Signature: *Dorinda Tiller* Date: 6/2/94

*184*

Food and Drug Administration  
Center for Devices and  
Radiological Health  
Office of Device Evaluation  
Document Mail Center (HFZ-401)  
1390 Piccard Drive  
Rockville, Maryland 20850

April 08, 1994

SANESE MEDICAL CORPORATION  
885 NORTHWEST BLVD.  
COLUMBUS, OH 43212  
ATTN: CHRISTOPHER N. SANESE

510(k) Number: K941507  
Received: 29-MAR-94  
Product: THEMO-FLO SYSTEM  
3  
LAPAROSCOPIC/HYS  
TEROSCOPIC

The Center for Devices and Radiological Health (CDRH), Office of Device Evaluation (ODE), has received the Premarket Notification you submitted in accordance with Section 510(k) of the Federal Food, Drug, and Cosmetic Act (Act) for the above referenced product. We have assigned your submission a unique 510(k) number that is cited above. Please refer prominently to this 510(k) number in any future correspondence that relates to this submission. We will notify you when the processing of your premarket notification has been completed or if any additional information is required.

The Safe Medical Devices Act of 1990 (SMDA), signed on November 28, states that you may not place this device into commercial distribution until you receive a letter from FDA allowing you to do so. Although the traditional timeframes for reviewing 510(k)s has been 90 days, it is now taking longer. These increasing response times have been caused by many factors, including a sharp increase in ODE's workload and increasingly complex device submissions. During 1992, we received about 1,500 more total submissions than we did the preceding year. We are troubled by these increases in response times and are making every effort to regain predictability in the timing of 510(k) reviews. Due to the increase in response times, CDRH has established a 510(k) Status Reporting System through which submitters may receive a status report on their 510(k) submissions(s) as follows:

- o Beginning 90 days after ODE receives your 510(k) submission, you may begin requesting status information. Submit requests via fax (301-443-8818) or via mail to:
  - 510(k) Status Coordinator
  - Division of Small Manufacturers Assistance (DSMA) (HFZ-220)
  - Center for Devices and Radiological Health, FDA
  - 5600 Fishers Lane
  - Rockville, Maryland 20857 USA

Because of staff limitations, we cannot answer telephone status requests.

- o 510(k) status requests should include:
  - (1) submitter's name and mailing address;
  - (2) requester's name, affiliation with the 510(k) submitter, mailing address, fax number (if applicable), telephone number, and signature; and

185

- (3) 510(k) information, including product name, 510(k) number, date logged in by ODE (as identified in acknowledgment letter from ODE), and name of contact person identified on firm's 510(k) submission.

Enclosed is a suggested format that you may use to ensure that you include all of the required information.

- o Within three working days after DSMA receives a submitter's status request, DSMA will send the submitter a fax or letter that includes:
  - (1) the branch to which the 510(k) has been assigned;
  - (2) the last action, and date of that action, that CDRH has taken regarding the 510(k), e.g., logging in an amendment, preparing a decision letter; and
  - (3) the position of the 510(k) in the reviewer's queue.

We request that 510(k) submitters make status inquiries no more than every four weeks. We do not have the resources to respond more frequently.

The SMDA also requires all persons submitting a premarket notification submission to include either (1) a summary of the safety and effectiveness information in the premarket notification submission upon which an equivalence determination could be based (510(k) summary), OR (2) a statement that safety and effectiveness information will be made available to interested persons upon request (510(k) statement). Safety and effectiveness information refers to information in the premarket notification submission, including adverse safety and effectiveness information, that is relevant to an assessment of substantial equivalence. The information could be descriptive information about the new and predicate device(s), or performance or clinical testing information. We cannot issue a final decision on your 510(k) unless you comply with this requirement.

Although FDA acknowledges that the law provides the 510(k) submitter an alternative, FDA encourages 510(k) submitters to provide a 510(k) statement to FDA and to make their safety and effectiveness information available to the public, excluding confidential manufacturing process information, in lieu of submitting a 510(k) summary to the agency until FDA promulgates a regulation on the content and format of 510(k) summaries. Since the law requires that FDA must make the 510(k) summary, or the source of information referred to in the 510(k) statement, publicly available within 30 days of making a substantial equivalence determination, we advise you that we may no longer honor any request for extended confidentiality under 21 CFR 807.95.

Additionally, the new legislation also requires any person who asserts that their device is substantially equivalent to a class III device to (1) certify that he or she has conducted a reasonable search of all information known, or otherwise available, about the generic type of device, AND (2) provide a summary description of the types of safety and effectiveness problems associated with the type of device and a citation to the literature, or other sources of information, upon which they have based the description (class III summary and certification). The

186

description should be sufficiently comprehensive to demonstrate that an applicant is fully aware of the types of problems to which the device is susceptible. If you have not provided this class III summary and certification in your premarket notification, please provide it as soon as possible. We cannot complete the review of your submission until you do so.

As of March 9, 1993, FDA has implemented the Good Manufacturing Practice(GMP) Pre-Clearance Inspection Program for all class III devices that are being reviewed under the premarket notification program. A letter of substantial equivalence cannot be sent until the finished device manufacturing site(s) and sterilization sites(s) as appropriate, have been identified and FDA has determined that the manufacturer(s) is in compliance with the GMP regulation (21 CFR Part 820).

Furthermore, the new legislation, section 522(a)(1), of the Act, states that if your device is a permanent implant the failure of which may cause death, you may be subject to required postmarket surveillance. If the premarket notification for your device was originally received on or after November 8, 1991, is subsequently found to be substantially equivalent to an Aneurysm Clip, Annuloplasty Ring, Artificial Embolization Device, Automatic Implanted Cardioverter Defibrillator System, Cardiovascular Intravascular Filter, Cardiovascular Permanent Pacemaker Electrode (Lead), Central Nervous System Fluid Shunt, Coronary Vascular Stent, Implantable Pacemaker Pulse Generator, Implanted Diaphragmatic/Phrenic Nerve Stimulator, Intracardiac Patch or Pledget, Intravascular Occluding Catheter, Replacement Heart Valve, Total Artificial Heart, Tracheal Prosthesis, Vascular Graft Prosthesis (less than 6 mm diameter), Vascular Graft Prosthesis (6 mm or greater diameter), Vena Cava Clip, or Ventricular Assist Device - Implant, you will be subject to the required postmarket surveillance and so notified of this determination in your substantially equivalent letter. (Some of the above listed types of devices may require a premarket approval application). This list is subject to change without notification. If you have any questions as to whether or not your device may be subject to postmarket surveillance or about this program, please contact the Postmarket Surveillance Studies Branch at (301) 594-0639.

Please note that the SMDA may have additional requirements affecting your device. You will be informed of these requirements as they become effective.

Please remember that all correspondence concerning your submission MUST be sent to the Document Mail Center (HFZ-401) at the above letterhead address. Correspondence sent to any address other than the Document Mail Center will not be considered as part of your official premarket notification submission. Because of equipment and personnel limitations we cannot accept telefaxed material as part of your official premarket notification submission, unless specifically requested of you by an FDA official.

187

If you have procedural or policy questions, please contact the Division of Small Manufacturers Assistance at (301) 443-6597 or their toll-free number (800) 638-2041, or contact me at (301) 594-1190.

Sincerely yours,

Marjorie Shulman  
Supervisory Consumer Safety Officer  
Premarket Notification Section  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

188

**PREMARKET NOTIFICATION (510(k)) STATUS REQUEST**

TO: 510(k) Status Coordinator  
Division of Small Manufacturers Assistance (HFZ-220)  
Center for Devices and Radiological Health, FDA  
5600 Fishers Lane  
Rockville, MD 20857  
USA  
Fax Number: (301) 443-8818

Please provide the status of the 510(k) identified below. Please send the information to the requester identified in section B by (check one):

fax  
 mail

**A. Sponsor Information:**

1. Name of 510(k) sponsor: \_\_\_\_\_

2. Sponsor's mailing address: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

**B. Requester information:**

1. Request name: \_\_\_\_\_

2. Requester affiliation with sponsor: \_\_\_\_\_

3. Requester mailing address: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

4. Request fax number (if applicable): \_\_\_\_\_

5. Requester telephone number: \_\_\_\_\_

**C. 510(k) information:**


1. Product name: \_\_\_\_\_

2. 510(k) number: \_\_\_\_\_

3. Date logged in by Office of Device Evaluation (ODE) (as identified in acknowledgment letter from ODE): \_\_\_\_\_

Name of contact person identified on firm's 510(k) submission: \_\_\_\_\_

.....  
I certify that the above information is accurate and truthful to the best of my knowledge.







**SANESE MEDICAL**  
CORPORATION

MANUFACTURER OF MEDICAL PRODUCTS

**K941507**

**RECEIVED**

29 MAR 94 13 06

FDA/CDRH/QDE/DMC

**OB**

**CLASS II DEVICE  
PREMARKET NOTIFICATION**

**CONFIDENTIAL - 807.95(b, c)  
CLAIMS MADE  
DO NOT LOG IN PUBLICLY  
AVAILABLE LOG**

March 18, 1994

Center for Devices and Radiological Health  
Document Mail Center (HFZ-401)  
Food and Drug Administration  
1390 Piccard Drive  
Rockville, MD 20850

**RE: 510(k) Notification of Sanese Medical Corporation's  
Intention to Market Its Thermo-Flo System 3  
Laparoscopic/Hysteroscopic Irrigation System and Sanese  
Disposable Irrigation/Aspiration Probe and Tubing.**

Dear Sir/Madam:

Pursuant to § 510(k) of the Food, Drug & Cosmetic Act, Sanese Medical Corporation (SMC) hereby submits premarket notification of its intention to market a laparoscopic/hysteroscopic irrigation pump system and disposable irrigation/aspiration probe and tubing. Enclosed is an original two copies of SMC's premarket notification for these devices.

Respectfully submitted,

Christopher N. Sanese  
Regulatory Correspondent  
Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, OH 43212  
(614) 291-8848

**190**

**SANESE MEDICAL**  
CORPORATION

MANUFACTURER OF MEDICAL PRODUCTS

PREMARKET NOTIFICATION  
OF CLASS II  
MEDICAL DEVICE

CONFIDENTIAL - 807.95(b, c)  
CLAIMS MADE  
DO NOT LOG IN PUBLICLY AVAILABLE  
LOG

March 16, 1994

Center for Devices and Radiological Health  
Document Mail Center (HFZ-401)  
Food and Drug Administration  
1390 Piccard Drive  
Rockville, MD 20850

RE: 510(k) Notification of Sanese Medical Corporation 's  
Intention to Market its Thermo-Flo System 3  
Laparoscopic/Hysteroscopic Irrigation System and  
Sanese Disposable Irrigation/Aspiration Probe and  
Tubing

Dear Sir/Madam:

I.

Premarket Notification of Intention to Market  
Laparoscopic/Hysteroscopic Irrigation System  
and Disposable Irrigation/Aspiration Probe  
and Tubing

Pursuant to 21 U.S.C. § 360(k) and 21 C.F.R. §§ 807.87 and 807.90, the Sanese Medical Corporation ("SMC") hereby submits premarket notification of its intention to market a Laparoscopic/Hysteroscopic Irrigation System and Disposable Irrigation/Aspiration Probe and Tubing, which will be marketed under the proprietary names Thermo-Flo System 3<sup>™</sup> and Sanese Disposable Irrigation/Aspiration Probe and Tubing.<sup>1/</sup> These are a Class II medical devices<sup>1/</sup> and are substantially equivalent in terms of safety and effectiveness to other laparoscopic irrigation pumps and aspiration probe and tubing devices currently on the market. (See Exhibits C and D.)

1 See 21 C.F.R. §§ 884.1690 and 1720. See also 21 C.F.R. § 880.5475 (Jet Lavage) and 21 C.F.R. § 878.4780 (Powered Suction Pump).

Center for Devices and Radiological Health  
March 16, 1994  
Page - 2 -

II.

The Section 510(k) Summary of Safety and Effectiveness  
Information for SMC's Thermo-Flo System 3 Laparoscopic/  
Hysteroscopic Irrigation System and Sanese Disposable  
Irrigation/Aspiration Probe and Tubing  
is attached as Exhibit A.

1992

III.

Supporting Material and Information

A. Names of Devices:

(i) Proprietary Names:

- (a) Thermo-Flo System 3 Laparoscopic/  
Hysteroscopic Irrigation System<sup>®</sup>
- (b) Sanese Disposable Irrigation/Aspiration Probe  
and Tubing™

(ii) Common, usual names:

- (a) Laparoscopic irrigation pump system
- (b) Disposable aspiration/irrigation probe and  
tubing.

(iii) Classification Names: Hysteroscope and  
Accessories and Gynecologic Laparoscope and  
Accessories, 21 C.F.R. §§ 880.1690 and 884.1720.

B. Regulatory Class:

Class II.

C. Actions Required for Compliance With Performance  
Standards: None. No standards have been  
promulgated to date.

D. Name and Address of Person Filing 510(k) Submission:

Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, OH 43212  
(614) 291-8848

E. Regulatory Correspondent of Person Filing 510(k)  
Submission:

Christopher N. Sanese  
Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, OH 43212  
(614) 291-8848

F. Registration Number of Device Manufacturer: The  
Company is not registered presently.

198

















A

102

SECTION 510(K) SUMMARY OF SAFETY AND  
EFFECTIVENESS INFORMATION FOR SMC'S  
THERMO-FLO SYSTEM 3 LAPAROSCOPIC/  
HYSTEROSCOPIC IRRIGATION SYSTEM AND  
SANESE DISPOSABLE IRRIGATION/ASPIRATION  
PROBE AND TUBING

Pursuant to Section 12(a) of the Safe Medical Devices Act of 1990, 21 U.S.C. § 360c(i)(3)A and 21 C.F.R. § 807.92, the following summary of information concerning the safety and effectiveness of SMC's Thermo-Flo System 3 Laparoscopic/Hysteroscopic Irrigation System and Sanese Disposable Irrigation/Aspiration Probe and Tubing is hereby submitted.

1. Name and address of device manufacturer submitting 510(k) notification:

Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, OH 43212  
(614) 291-8848

2. Regulatory correspondent.

Christopher N. Sanese  
Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, OH 43212  
(614) 291-8848

3. Date summary was prepared: March 17, 1994

4. Names of Devices:

(i) Proprietary Names:

- (a) Thermo-Flo System 3 Laparoscopic/  
Hysteroscopic Irrigation System<sup>®</sup>
- (b) Sanese Disposable Irrigation/Aspiration Probe  
and Tubing<sup>™</sup>

(ii) Common, usual names:

- (a) Laparoscopic irrigation pump system
- (b) Disposable aspiration/irrigation probe and  
tubing.

(iii) Classification Names: Hysteroscope and  
Accessories and Gynecologic Laparoscope and  
Accessories, 21 C.F.R. §§ 880.1690 and 884.1720.

5. Predicate devices to which SMC is claiming substantial equivalence:

202

SMC claims substantial equivalence for purposes of section 510(k) only to:

A. Irrigation Pump System:

The Niagra™ High-Flow Irrigator System, 510(k) Number K924530 by Cabot Medical, 2021 Cabot Boulevard, West Langhorne, PA 19047 (800-523-6078).

B. Disposable Probe and Tubing:

(1) Probe:

Corson™ Disposable Suction Irrigation Probe by Cabot Medical, 2021 Cabot Boulevard, West Langhorne, PA 19047 (800-523-6078).

(2) Tubing:

Core Dynamics Inc. Disposable Irrigation/Suction Tubing Set (510(k) Number K910246) by Core Dynamics, Inc., 9951 Atlantic Boulevard, Jacksonville, Florida 3225 (904-727-0910).

6. Description of the devices and primary principles of operation:

A. Description of irrigation pump system.

This device consists of a fabric enclosure (pressure envelope) in which a bag of irrigation fluid is placed. The fabric enclosure contains an inflatable bladder that is pneumatically inflated to pressurize the fluid in the irrigation fluid bag. Inflation of the bladder is provided by an air compressor which has a maximum voltage leakage current of 100 microamps. This complies with Underwriters Laboratory Standard 544.

The pressure envelope also contains a heating element to assist in maintaining the temperature of preheated irrigation solution while being used in the pressure envelope. There is also a thermal regulator to control the temperature of the heating element. The heating element and the thermal regulator are powered by an Underwriter Laboratory approved 24 volt electrical power supply.

The air compressor and heating element are controlled by a microcontroller. These components are controlled by the micro-controller through activating solid state switches which control electrical power to the pump, pump solenoids and the heater.

203

The pressurized preheated fluid is then introduced into the endoscopic surgery site through a medical grade plastic tube and probe that are discussed below.

B. Description of disposable probe and tubing:

The probe consists of a medical grade plastic, activation valve and metal shaft that will deliver and aspirate fluids to and from an operative site.

The tubing consists of two 15 foot length medical grade surgical tubing that are approximately 10 millimeters in diameter which would be connected to the probe and the irrigation solution bag and the pump when used for irrigation or to the probe and the suction unit when used for aspiration. The probe and tubing are made of medical grade plastic and metal that have been determined by FDA to be safe and effective for this intended usage.

C. Primary Principles of Operation:

(i) Irrigation Pump System:

The primary principles of operation are pressurization through an electrically powered pneumatic compressor for pressurized irrigation fluid flow and thermal transfer through an electrically-powered thermo-conducted heat transfer unit to maintain preheated fluid.

(ii) Probe and Tubing:

With respect to the aspiration/irrigation use of the probe and tubing, the principle operation is vacuum and pressurization.

7. Statement of the purpose for which the device will be recommended:

These devices are intended to be used for irrigation and aspiration during adult and pediatric hysteroscopic or laparoscopic (endoscopic) surgery.

8. Biocompatibility and effectiveness of the device:


(i) Biocompatibility:

The only aspect of the device that comes in contact with the patient is the distal end of the probe. The probe and the probe tubing have been subject to biocompatibility tests and are safe and effective for such usage. The irrigation pump has a maximum voltage

204

leakage current of 100 microamps in accordance with U.L. Standard 544. Performance tests of the pump demonstrate that it is safe and effective for its intended use. Also, all components of the pump meet U.L. Standards.

Respectfully submitted,



Christopher N. Sanese  
Regulatory Correspondent  
Sanese Medical Corporation  
885 Northwest Boulevard  
Columbus, OH 43212  
(614) 291-8848

705




B

2060

IRRIGATION PUMP SYSTEM

**THERMO-FLO**  
SYSTEM 3



SANESE MEDICAL  
CORPORATION

The new Thermo-Flo System 3 was developed to eliminate the major problems faced everyday by the surgeon and o.r. unit manager.

207

The Thermo-Flow System 3 was developed to eliminate the major problems faced everyday by the surgeons and o.r. unit managers. These are: Cross contamination, inadvertent patient hypothermia, inadequate irrigation solution capacity, ease of set up, tube sets and disposable probes that leak, noise and physical space. Sanese Medical answers these problems with the Thermo-Flow System 3.

1) **CROSS CONTAMINATION** occurs when non sterile fluid enters the patients abdominal cavity through the irrigation line. The Thermo-Flow System 3 eliminates cross contamination.

2) **INADVERTENT PATIENT HYPOTHERMIA** is the result of the combination of gas and fluids ranging in temperature from 65 to 72 degrees farenhite being introduced into the patient. The Thermo-Flow System 3 maintains the solution temperature at 100 degrees farenhite. Though this system alone cannot prevent hypothermia, it will dramatically effect the reduction of the core temperature due to the existing conditions of surgery.


3) **SOLUTION CAPACITY / EASE OF SET UP**  
The Thermo-Flow System 3 meets the three liter needs of todays and tomorrows laparoscopic procedureds. This amount provides less disruption of the procedure, by not replacing the solution, cap and tube inserts that promotes the inadvertent contamination by the person replacing the solution as often.

**THERMO-FLOW SYSTEM 3 PRODUCT SYSTEM**

- 1) The Thermo- Flow System 3 Irrigation Pump.....TFS3 - one per case
- 2) SMC / TFS-3 Tube set with disposable probe..... TFS3 - 1-WP one per case  
(15'tube set: aspiration included)
- 3) SMC / TFS-3 Tube set without disposable probe..... TFS3 - 1-WO ten sets per case  
(15' tube set: aspiration included)
- 4) SMC-NS Tube set with disposable probe.....TFS3 - N - 1WP ten per case  
(15' tube set: aspiration included)
- 5) SMC-NS Tube set without disposable probe..... TFS3 - N -1W0 ten sets per case  
(15' Tube set: aspiration included)

**KEY ADVANTAGES OF THE SMC / TFS3**

- \*Delivers strong stream of irrigation fluid
- \*Choice of disposable or reusable probe
- \*No capitol equipment necessary
- \*Simplified O.R. set up / uses I.V.
- \*Avoids equipment purchase
- \*Low level sound
- \* 24 hour, 7 days a week service available



**SANESE MEDICAL**  
CORPORATION

885 Northwest Blvd. / Columbus, Ohio 43212

For further information or to order call (614) 291-8848

Fax (614) 291-1389

LOT NO. 0300



**CONTENTS:**

1 - MODEL TFS 3 LAPAROSCOPIC / HYSTEROSCOPIC  
IRRIGATION SYSTEM

**INDICATED FOR IRRIGATION DURING LAPAROSCOPIC &  
HYSTEROSCOPIC SURGERY**

**Distributed by:**

**SANESE MEDICAL**  
CORPORATION

885 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-8848 FAX (614) 291-1389

**CAUTION: Federal ( U.S.A. ) Law restricts this device to the sale by,  
or on the order of a physician.**

*Non - Sterile - Recommended for irrigation in adult and pediatric  
laparoscopic / hysteroscopic Surgery*

209

**CONTENTS: MODEL NO. TFS3 - N 2WP**

1 - Single use, disposable 15' irrigation / aspiration tube set with probe

**LOT NO. 0500**

**FOR SINGLE USE ONLY**

**CONTENTS STERILE IN UNOPENED, UNDAMAGED PACKAGES.**

**INDICATED FOR IRRIGATION / ASPIRATION DURING LAPAROSCOPIC  
HYSTEROSCOPIC SURGERY.**

***Distributed by:* SANESE MEDICAL**  
CORPORATION

885 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-8848 FAX (614) 291-1389

*RECOMMEND FOR ADULT AND PEDIATRIC  
LAPAROSCOPIC / HYSTEROSCOPIC SURGERY*

*CAUTION: Federal ( U.S.A.) law restricts this device to sale by, or on the order of a physician.*

710

**CONTENTS: MODEL NO. TFS3 - 3 WOP**

1 - Single use, disposable 15' irrigation / aspiration tube set.

**LOT NO. 0500**

**FOR SINGLE USE ONLY**

**CONTENTS STERILE IN UNOPENED, UNDAMAGED PACKAGES.**

**INDICATED FOR IRRIGATION / ASPIRATION DURING LAPAROSCOPIC  
HYSTEROSCOPIC SURGERY.**

***Distributed by:*** **SANESE MEDICAL**  
CORPORATION

885 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-8848 FAX (614) 291-1389

*RECOMMEND FOR ADULT AND PEDIATRIC  
LAPAROSCOPIC / HYSTEROSCOPIC SURGERY*

CAUTION: *Federal ( U.S.A.) law restricts this device to sale by, or on the order of a physician.*

211

**CONTENTS: MODEL NO. TFS3 - 3 WP**

1 - Single use, disposable 15' irrigation / aspiration tube set with probe

**LOT NO. 0500**

**FOR SINGLE USE ONLY**

**CONTENTS STERILE IN UNOPENED, UNDAMAGED PACKAGES.**

**INDICATED FOR IRRIGATION / ASPIRATION DURING LAPAROSCOPIC  
HYSTEROSCOPIC SURGERY.**

***Distributed by:*** **SANESE MEDICAL**  
CORPORATION

885 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-8848 FAX (614) 291-1389

*RECOMMEND FOR ADULT AND PEDIATRIC  
LAPAROSCOPIC / HYSTEROSCOPIC SURGERY*

CAUTION: Federal ( U.S.A.) law restricts this device to sale by, or on the order of a physician.

212

**OPERATIONS MANUAL**



**LAPAROSCOPIC / HYSTEROSCOPIC IRRIGATION  
SYSTEM**

**SANESE MEDICAL**  
CORPORATION

213

885 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-8848 FAX (614) 291-1389



## ***SAFETY PRECAUTIONS***

**\*\* USE ONLY WITH TFS 3 IRRIGATION TUBE SETS \*\***

- \* Do not use in volatile atmosphere.
- \* Use only with 120 volt grounded receptacle.
- \* Do not immerse unit.

## **READ BEFORE ADJUSTING AND OPERATING UNIT**

It is recommended that all O.R. Unit Managers and others responsible for the operation of this unit become thoroughly familiar with its capabilities and proper operating procedures prior to actual patient use.

## ***PRODUCT DESCRIPTION***

The **THERMO-FLO SYSTEM 3<sup>TM</sup>**\* LAPAROSCOPIC / HYSTEROSCOPIC IRRIGATION SYSTEM consists of the following components:

- 1 - **THERMO - FLO SYSTEM 3<sup>TM</sup>** control device
- 1 - **THERMO - FLO SYSTEM 3<sup>TM</sup>** pressure enclosure device
- 2 - Transportation pouch

\*U. S. and International patent pending  
**SANESE MEDICAL**  
CORPORATION

885 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-8848 FAX (614) 291-1389

214

## ***SPECIFICATIONS***

### **MODEL TFS3\***

#### **POWER**

120 V AC 1 AMP OR 220 V AC 0.5 AMP 50/60 Hz  
( 95 to 250 V AC 47 to 63 Hz

#### **WEIGHT**

7 LB. ( with out stand )

#### **PRESSURE RANGE**

100 mm TO 800 mm HG

#### **HEATER TEMPERATURE RANGE**

80 F to 120 F  $\pm$  2 F

Safety Shutdown Temperature: 125 F

#### **MAXIMUM LEAKAGE CURRENT**

Control unit: 100 micro amps ( per UL 544 )

Pressure Enclosure: 10 micro amps (per UL544 )

Maximum Ground Impedance: 0.1 Ohm ( per UL 544 )

Fuse. Slow Blow 2.0 Amp. ( Located on back of control unit )

#### **ENVIRONMENTAL OPERATING CONDITIONS**

40 F TO 150 F, Humidity: non-condensing

\* U. S. and International patent pending

Specifications subject to change without notice.

215

**SANESE MEDICAL**  
CORPORATION

## ***SAFETY FEATURES***

The **THERMO-FLO SYSTEM 3** IRRIGATION SYSTEM was designed with the following safety features:

### **OVER PRESSURE**

In the event of an over pressure condition a pressure switch cuts power to the air supply.

### **OVER TEMPERATURE**

In the event of an over temperature condition, an independent software circuit cuts the power supply to the heater.

### **LEAKAGE CURRENT**

The leakage current of the pressure cuff is less than 10 micro amps.

### **WATCHDOG TIMER**

If the micro controller is not executing code properly, the watchdog timer resets the unit to its initial conditions.

**SANESE MEDICAL**  
CORPORATION

885 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-8848 FAX (614) 291-1389

## ***OPERATING THE TFS3™ IRRIGATION SYSTEM***

### *DESCRIPTION OF CONTROLS*

**Power Switch:** The power switch is a double pole single throw 10 Amp switch. (Both the line and neutral are switched per UL544).

**Pressure Adjust Switch:** Located to the left and below the display is a three position rocker switch which is used to adjust the pressure. The center position is off. The top (Up arrow) and bottom (Down arrow) positions are momentary. Pressing the Up arrow increases the fluid pressure. Pressing the Down arrow decreases the fluid pressure.

**Temperature Adjust Switch:** Located to the right and below the display is a three position rocker switch which is used to adjust the heater temperature. The center position is off. The top (Up arrow) and bottom (Down arrow) positions are momentary. Pressing the Up arrow increases the heater temperature. Pressing the Down arrow decreases the heater temperature.

**Load / Run Switch:** This is a two position rocker switch. When in the LOAD position the pump is turned off and the solenoid valve is opened allowing an empty bag of fluid to be unloaded and a full bag loaded. When in the RUN position the pump and solenoid are under control of the microcontroller which turns on the pump and opens the solenoid as needed to keep the pressure at the setpoint.

**Display:** The display is a 16 character single line LCD Module operated by the microcontroller. In normal operation it displays the pressure and temperature setpoints. It is also used to prompt the user.

**Fuse:** A 2 amp fuse is located on the back of the unit.

**Air Connector:** On the back of the unit is located a Colder Connectors into which plugs the mating connector located on the end of the air bladder tube.

**Heater Electrical Connector:** On the back of the controller is located a connector used to electrically connect the heaters and temperature sensors to the controller.

**SANESE MEDICAL**  
CORPORATION

## *DESCRIPTION OF OPERATION*

The controller is plugged into a standard wall outlet. The power switch is turned on and the LOAD/RUN switch is placed in the LOAD position. A fresh bag of irrigation fluid is placed in the pressure enclosure. The irrigation bag is connected to the tubing set / probe through the stick fittings which protrude through a hole in the bottom of the pressure enclosure. The LOAD/RUN switch is placed in the RUN position. When power is first turned on the pressure and temperature settings are at minimum. The user adjusts the Pressure and Temperature setpoints to their desired settings. The user then depresses the trumpet valve on the probe to obtain irrigation fluid. After a time the bag becomes empty. The LOAD/RUN switch is then moved to LOAD and the empty bag is replaced with a full one. Then the LOAD/RUN switch is moved to RUN. Because the machine has not been turned off the pressure and temperature settings remain at the value they had before the new bag of fluid was added.

812

**SANESE MEDICAL**  
CORPORATION

885 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-8848 FAX (614) 291-1389

USAGE RECOMMENDATION

The **THERMO-FLO SYSTEM 3™** was designed for Laparoscopic / Hysteroscopic surgery and is only to be used by licensed, certified surgical personnel.

RECOMMENDED CARE

Do not leave unit plugged in.

Do not use abrasive cleaning and disinfecting solutions.

LIMITED WARRANTY

Sanese Medical Corporation ("Sanese") warrants the **THERMO-FLO SYSTEM 3™** to be free from defects in material and workmanship for a period of one (1) year from date of initial purchase. No warranty shall apply if the goods have been damaged by accident, abuse, misuse, or misapplication, or as a result of service or modification by other than a person authorized by Sanese. Sanese's liability under this warranty, and buyer's exclusive remedy is limited to the cost of materials and labor to repair defective goods, or to their replacement, and in no event shall exceed the purchase price. Repair or replacement of defective goods under warranty will be made only upon their return to Sanese after notice to Sanese and Buyer's receipt of shipping instructions. Upon receipt of goods returned under warranty, Sanese will inspect the goods, will notify buyer of the extent of repair or replacement which Sanese will perform under warranty, which shall be conclusive of Sanese's liability.

Sanese is not responsible for incidental consequential damages resulting from the breach of any expressed or implied warranty, including damage to property and to the extent permitted by law, damages for personal injury. The warranty contained herein is in lieu of all warranties, expressed or implied. No statement of any representative shall extend Sanese's liability as herein established or limited.

**SANESE MEDICAL**  
CORPORATION

885 NORTHWEST BLVD. COLUMBUS, OHIO 43212 (614) 291-8848 FAX (614) 291-1389

219  
612

C

220

# Niagara™ High-Flow Irrigator Systems

For quiet, convenient, cost-effective irrigation, the Niagara irrigator sets the standard.


Niagara Irrigator's exclusive pressurized-chamber technology delivers flow rates up to 3L/minute—nearly 50% higher than disposable diaphragm pumps—for confident visualization and precise hydrodissection.

Niagara pumps are easy to use. Just hang from any IV pole...attach to OR air, nitrogen, or O<sub>2</sub>, or bottled CO<sub>2</sub>...hang and spike the solution bags...close the chamber...and you're ready.

Niagara Irrigators are built for safety—with double-layer protection to ensure that the pressurizing gas never contacts the irrigant.


And, because Niagara Irrigators are reusable procedure after procedure, they cut pump/probe costs nearly in half.

Choose from three versatile models designed to fit a wide range of needs.




## 1L Single-Pump Niagara Irrigator

For low-volume applications, the Niagara 1L Single-Pump model provides high-flow performance in a compact, cost-effective package.



## 2L Dual-Pump Niagara Irrigator

Equally at home in long and short procedures, the Niagara 2L Dual-Pump accommodates two separately pressurized 1L fluid bags—with instant switchover to assure uninterrupted availability of irrigant.




NIAGARA™

3L High Flow Irrigator

CAUTION

- Do Not Use for Infusion of Blood or IV's
- Do Not Use for Hysteroscopy
- Do Not Exceed 500mm Hg Pre-stare







Air/Nitrogen Connector




Oxygen Connector



CO<sub>2</sub> Tank Connector



Irrigation Pump Stand




**Cabot Irrigation Pump**

For use with bottled solutions, the CO<sub>2</sub>-powered Cabot Irrigation Pump offers continuously adjustable pressure for a range of flow rates—up to 2L/minute at pressures up to 775mm Hg—and continuous availability of irrigant, thanks to its two-bottle design.


(See last page for complete selector)

*Nobody puts it together  
like Cabot.*

**800-523-6078**




*Tissue Collector*



cabot medical

# LAPAROSCOPIC INSTRUMENTS

## Corson™ Disposable Suction/ Irrigation Probe without Tubing for Use with Disposable Pumps




005407-901

- Can position suction and irrigation functions according to physician preference.
- Instrument channel provides ability to dissect, lase, cut, cauterize, irrigate or aspirate through a single port.
- High flow design allows for maximum irrigation and suction including aspiration of blood clots and gallstones.
- Innovative trumpet valve virtually eliminates clogging, leaking and sticking.

## CORSON™ DISPOSABLE SUCTION/ IRRIGATION PROBE WITHOUT INTEGRATED SUCTION AND IRRIGATION TUBING

(Package of 10)  
28 cm

005407-901



005764-902

## Corson™ Probe Retractors

The Corson™ Probe Retractors are modified versions of the Corson™ Probes and are available in 28 cm working lengths. They are manufactured without holes in the distal ends to enable the probes to be used to suction grasp/retract small organs.

## CORSON™ PROBE RETRACTORS

(Package of 8)

Retractor

005765-902

Retractor with Spike

005764-902

## Big Sucker

10mm Tissue Aspiration and Collection System

Aspirate and collect stones, large clots, biopsies, and ectopic pregnancies.

- Large bore, straight-through design permits aspiration of both rigid and soft specimens.
- Suction-relief collar allows quick release of oversize tissues.
- Tissue collector preserves intact specimens for histopathology.
- Can be used simultaneously with 5mm suction/irrigation systems.



005773-901

Big Sucker

005773-901

# The Big Sucker 10mm Tissue Aspiration & Collection System


This remarkable new system extends the range of laparoscopic surgery with unprecedented control.

Now you can easily and safely collect "impossible" tissues—including large clots, multiple stones, biopsies, and ectopic pregnancies—without time-consuming, tissue-destroying dissection.


The large-bore cannula and straight-through valve design permit rapid, high-suction aspiration of both rigid and friable specimens while the exclusive tissue collector captures and preserves specimens intact for histopathology.

The Big Sucker System maximizes safety and control—with a suction-relief collar that guards against unwanted tissue aspiration, allows quick release of oversize tissues and a clear cannula that lets you see exactly what is being aspirated.


And it maximizes convenience by permitting simultaneous use of 5mm suction/irrigation systems.




*Y-Connector for 5mm Suction/  
Irrigation Hook-up*



*Suction-Relief Collar*



*Straight-Through Valve Design*



*Large-Bore Transparent Cannula*

2

226



## DISPOSABLE IRRIGATION/SUCTION TUBING SET

For Laparoscopy, Pelviscopy and Cholecystectomy procedures

### Performance Features

Fits all irrigation/probes with trumpet valves \*

Easy to attach luer lock connectors to prevent fluid or vacuum leakage \*

Universal spike on irrigation tubing attaches to any solution bag \*

Pinch clamp on irrigation tubing allows leak proof disconnection \*


Light weight and flexible for better probe control \*

Non-kinking \*

8 Foot lengths increases accessibility \*

Sterile package \*


One time use \*



*A cost effective disposable irrigation/suction tubing set for all endoscopic procedures. When used directly with a trumpet valve irrigation/suction probe, costly irrigation pumps are eliminated. Easy to set up and minimal equipment monitoring is needed.*

289

Follow diagram below for usage and assembly



distributor:

**Vector Medical Corporation**  
 5859 South Garnett  
 Tulsa, OK 74146  
 (918) 252-5400

**ORDER INFORMATION**

ORDER #	DESCRIPTION	PACKAGE
<b>CD3400A</b>	<i>Disposable Irrigation / Suction tubing</i>	<b>6 per box</b>

**CORE DYNAMICS,™ INC.**

9951 Atlantic Boulevard Jacksonville, Florida 32225 (904)727-0910 fax (904)724-7077

228

© CORE Dynamics, Inc. 1991

manufactured in the U.S.A

patent pending


printed in the U.S.A 3/91

# LAPAROSCOPIC INSTRUMENTS

## Corson™ Disposable Suction/Irrigation Probe for Bottle Fluids

This probe is specifically designed to be compatible with a broad range of Cabot Medical instrument inserts that dissect, lase, cut, cauterize, irrigate or aspirate.

- High-flow design allows for maximum irrigation and suction, including aspiration of blood clots and gallstones.
- Innovative trumpet design virtually eliminates clogging, leaking and sticking.



004946-901

## CORSON™ DISPOSABLE SUCTION/IRRIGATION PROBE WITH INTEGRATED SUCTION AND IRRIGATION TUBING


*(Package of 8 for Bottle Fluids)*

24 cm	004946-901
28 cm	004946-902
28 cm without holes	005765-902
32 cm	004946-903

## Corson™ Disposable Suction/Irrigation Probe with Spike for Bag Fluids

This probe is an adaptation of the Corson™ Probe for procedures using the Niagara™ Irrigator or gravity flow for irrigation.

- Dual spike allows irrigant to be delivered via multi-bag setup.
- Customized check valves prevent fluid from flowing backward into irrigation bags.



005200-901

## CORSON™ DISPOSABLE SUCTION/IRRIGATION PROBE WITH SPIKE WITH INTEGRATED SUCTION AND IRRIGATION TUBING

*(Package of 8 for Bag Fluids)*


24 cm	005200-901
28 cm	005200-902
28 cm without holes	005764-902
32 cm	005200-903



# LAPAROSCOPIC INSTRUMENTS

## Corson™ Disposable Suction/ Irrigation Probe without Tubing for Use with Disposable Pumps

- Can position suction and irrigation functions according to physician preference.
- Instrument channel provides ability to dissect, lase, cut, cauterize, irrigate or aspirate through a single port.
- High flow design allows for maximum irrigation and suction including aspiration of blood clots and gallstones.
- Innovative trumpet valve virtually eliminates clogging, leaking and sticking.




005407-901

## CORSON™ DISPOSABLE SUCTION/ IRRIGATION PROBE WITHOUT INTEGRATED SUCTION AND IRRIGATION TUBING

(Package of 10)  
28 cm

005407-901



005764-902

## Corson™ Probe Retractors

The Corson™ Probe Retractors are modified versions of the Corson™ Probes and are available in 28 cm working lengths. They are manufactured without holes in the distal ends to enable the probes to be used to suction/grasp/retract small organs.

## CORSON™ PROBE RETRACTORS

(Package of 8)

Retractor

005765-902

Retractor with Spike

005764-902

## Big Sucker


10mm Tissue Aspiration and Collection System

Aspirate and collect stones, large clots, biopsies, and ectopic pregnancies.

- Large bore, straight-through design permits aspiration of both rigid and soft specimens.
- Suction-relief collar allows quick release of oversize tissues.
- Tissue collector preserves intact specimens for histopathology.
- Can be used simultaneously with 5mm suction/irrigation systems.

Big Sucker

005773-901



005773-901

25

251

December 9th., 1993

**Mr. Christopher N. Sanese**  
**President**  
**Sanese Medical Corporation**  
**885 Northwest Blvd.**  
**Columbus, Ohio 43212**

**Dear Chris,**

**This letter is a follow-up to our conversation on Monday, December 6th. The USP Class VI information that you requested is not readily available from the various vendors that will be providing materials for the tubing set. The normal response to my inquiry was "We do not perform routine Class VI testing to the materials used in the production of our products." As a matter of fact, we at Surgimedics can make the same claim. Class VI testing is time consuming and expensive and is not always a requirement of our customers. I have included a test report that I did obtain from Natvar for tubing that could possibly be used in the tubing set.**

**Following is a listing of parts, materials and vendors for the components that I have selected for use in the tubing set.**

**Tubing: Clear PVC (70 durometer Shore "A") Purchased from Natvar Co. Formulation #660.**

**Male Luer: Natural Acrylic CYRO Industries G-20 Hi-Flo. Purchased from Kippmed #1001-187-004.**

**Spike Cover: Low Density Polyethylene, Eastman, Tenite 1830 "A" or "E" (Natural)  
Purchased From McGaw #X52-018.**

**Suction Connector: Blue PVC. Purchased from Qosina #57015**

**Non-Vented Spike: The material that I selected for the custom molded Non-Vented Spike is a White Polycarbonate GE Lexan #144R-112. We currently use this material to produce various medical devices. Attached is a copy of a letter that I received from General Electric concerning regulatory issues.**

**Chris, I hope that this is enough information to get you started with your 510(k) filing. I will continue to explore the possibility of obtaining Class VI information on the above items but I expect very limited success.**

**cc: H. Dutra**  
**C. Meehan**

Surgimedics/TMP  
9287 Smucker Road  
Orrville, Ohio 44667  
216/669-2271

**Sincerely,**

  
**Roger D. Ellis**  
**Manufacturing Engineer**





















# Code of federal regulations

## Food and Drugs

# 21

PARTS 170 TO 199  
Revised as of April 1, 1988



### § 177.1580 Polycarbonate resins.

Polycarbonate resins may be safely used as articles or components of articles intended for use in producing, manufacturing, packing, processing, preparing, treating, packaging, transporting, or holding food, in accordance with the following prescribed conditions:

(a) Polycarbonate resins are polyesters produced by:

(1) The condensation of 4,4'-isopropylidenediphenol and carbonyl chloride to which may have been added certain optional adjuvant substances required in the production of the resins; or by

(2) The reaction of molten 4,4'-isopropylidenediphenol with molten diphenyl carbonate in the presence of the disodium salt of 4,4'-isopropylidenediphenol.

(3) The condensation of 4,4'-isopropylidenediphenol, carbonyl chloride, and 0.5 percent weight maximum of  $\alpha, \alpha$ -bis (6-hydroxy-*m*-tolyl) mesitol to which may have been added certain optional adjuvant substances required in the production of branched polycarbonate resins.

(b) The optional adjuvant substances required in the production of resins produced by the methods described in paragraph (a)(1) and (3) of this section may include substances generally recognized as safe in food, substances used in accordance with a prior sanction or approval, and the following:

List of substances	Limitations
<i>p</i> -tert-Butylphenol .....	
Chloroform .....	
Ethylene dichloride .....	
Heptane .....	
Methylene chloride .....	
Monochlorobenzene .....	Not to exceed 500 p.p.m. as residual solvent in finished resin.
Phenol (CAS Reg. No. 106-95-2) .....	
Pyridine .....	
Toluene: (CAS Reg. No. 106-88-3) .....	Not to exceed 800 parts per million as residual solvent in finished resin.
Triethylamine .....	

(c) Polycarbonate resins shall conform to the specification prescribed in paragraph (c)(1) of this section and shall meet the extractives limitations prescribed in paragraph (c)(2) of this section.

(1) *Specification.* Polycarbonate resins can be identified by their characteristic infrared spectrum.

(2) *Extractives limitations.* The polycarbonate resins to be tested shall be ground or cut into small particles that will pass through a U.S. standard sieve No. 6 and that will be held on a U.S. standard sieve No. 10.

(i) Polycarbonate resins, when extracted with distilled water at reflux temperature for 6 hours, shall yield total extractives not to exceed 0.15 percent by weight of the resins.

(ii) Polycarbonate resins, when extracted with 50 percent (by volume) ethyl alcohol in distilled water at reflux temperature for 6 hours, shall yield total extractives not to exceed 0.15 percent by weight of the resins.

(iii) Polycarbonate resins, when extracted with *n*-heptane at reflux temperature for 6 hours, shall yield total extractives not to exceed 0.15 percent by weight of the resins.

[42 FR 14572, Mar. 15, 1977, as amended at 46 FR 23227, Apr. 24, 1981; 49 FR 4372, Feb. 6, 1984; 50 FR 14096, Apr. 10, 1985]

ZUR

F

243

(b)(4) Trade Secret Process - Testing



**Objective:**

To determine the pressure at which the pressure enclosure fails in order to establish that the safety factor is sufficient.

**Conditions:**

(b)(4) Trade Secret Process - Testing



**Results:**

(b)(4) Trade Secret Process - Testing



**Discussion:**

(b)(4) Trade Secret Process - Testing



2004

(b)(4) Trade Secret Process - Testing

**Objective:**

To run the irrigation pump and pressure enclosure at rated pressure for (b) and observe performance.

**Conditions:**

(b)(4) Trade Secret Process - Testing

**Results:**


(b)(4) Trade Secret Process - Testing

SH2




G

246




U2




248

# CONTROLLER / FABRIC ENCLOSURE CONNECTIONS LABEL (LOCATED ON BACK OF CONTROLLER)




249

# CONTROLLER RATINGS / WARNINGS LABEL



250



<b>A</b>	-----	<b>AS RELEASED</b>	<b>9/12/93</b>	<b>KAK</b>
REV	ECN No.	REVISION	DATE	DWN
<b>MATERIAL:</b>	-----	<b>NEXT LEVEL:</b>	<b>Sanese Medical</b>	
<b>FINISH:</b>	-----	<b>MODEL:</b>		
<b>UNLESS OTHERWISE SPECIFIED ALL THE DIMENSIONS ARE IN INCHES.</b>			<b>TITLE:</b>	
			<b>Electrical Block Diagram</b>	
<b>DWN: KEVIN KELLY</b>		<b>SCALE:</b>	<b>SHT</b>	<b>OF</b>
<b>DATE: 9/12/93</b>		<b>NONE</b>	<b>1</b>	<b>1</b>
			<b>DWG/PART NO.</b>	<b>SIZE</b>
			<b>Elec_1</b>	<b>A</b>

251






























































































H


791



# SANESE PROBE



*2912*



793

(b)(4) Trade  
Secret Process -  
Engineering

(b)(4) Trade  
Secret Process  
- Engineering

(b)(4) Trade  
S t (b)(4) Trade  
Secret Proc (b)(4)  
Trade  
S t

(b)(4) Trade  
Secret  
P

(b)(4)  
Trade  
Secret

**OVERALL DIMENSIONS**

**PROBE & TUBE SET STERILIZATION METHOD**

\* **Gamma Radiation:** as follows;

**STERILITY ASSURANCE LEVEL - S.A.L. =  $10^{-6}$**

**Dose to be determined by AAMI Method 1,**

**Dose to be not less than 2.5 M Rads**

**Materials**

**Materials for the Probe**

**Probe Plastic Valve Housing**

GE Lexan -GR or Cyloloc ABS HP or Equivalent

Use Methylene Chloride as the bonding agent

**Probe Tube**

302 Stainless Steel Seamless Tube

**Tubing and Pressure Seals**

Norton Tygon S-50-HL PVC Shore A,64