

MAY - 6 1997

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Mr. Jeffrey R. Mannion
Regulatory Affairs Associate
AngioDynamics, Inc.
603 Queensbury Avenue
Queensbury, New York 12804

Re: K964034

Transjugular Access Set Regulatory Class: II (two)

Product Code: DYB

Dated: February 4, 1997 Received: February 5, 1997

Dear Mr. Mannion:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the <u>Code of Federal Regulations</u>, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in

regulatory action. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for <u>in vitro</u> diagnostic devices), please contact the Office of Compliance at (301) 594-4648. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597.

Sincerely yours,

Thomas J. Callahan, Ph.D.

Director

Division of Cardiovascular, Respiratory,

and Neurological Devices
Office of Device Evaluation

Thomas J. Cellelon

Center for Devices and

Radiological Health

Enclosure

Page <u>1 of 1</u>

510(k) Number (if known): <u>K964034</u>
Device Name: <u>Transjugular Access Set</u>
Indications For Use:
The AngioDynamics, Inc. Transjugular Liver Access Set is intended to be used for transjugular liver access during diagnostic and interventional procedures.
ivision Sig vision of ? Neurolog  Oro(k) Number
vision Sign-Off) vision of Cardiovascular, Respiratory, Neurological Devices  (k) Number
(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)  Concurrence of CDRH, Office of Device Evaluation (ODE)
Prescription Use or

### 510(K) ROUTE SLIP


510/k) NUMBER	K964034 PANEL CV	DIVISION DCRND BRANCH ICDG
		Γ
	CATHETER INTRODUCER	
		r D
PRODUCT CODE	DYB INTRODUCER, CAIREL	ER
SHORT NAME CONTACT DIVISION ADDRESS PHONE NO.	ANGIODYNAMICS BRIAN KUNST  603 QUEENSBURY AVE. QUEENSBURY, NY 12804 (518) 798-1215	FAX NO. ( <u>518</u> ) <u>798-3625</u> REGISTRATION NO. <u>1319211</u> 1222273
DATE RECEIVED	MISSION <u>30-SEP-96</u> IN ODE <u>08-OCT-96</u> ECISION	DATE DUE TO 510(K) STAFF 22-DEC-96  DATE DECISION DUE 06-JAN-97  DECISION DATE
SUPPLEMENTS		DUE POS DUE OUT
<u>S001</u>	<u>04-FEB-97</u> <u>05-FEB-97</u>	<u>21-APR-97</u> <u>06-MAY-97</u>
CORRESPONDENCI	E SENT DUE BACK	
C001	06-JAN-97 05-FEB-97	HOLD LETTER
	) identified as a Class	
		MAY 6 1997

Records processed under FOIA Request #2018-3544; Released by CDRH on 09 13-2018

		•			Memorandum
	From:	Reviewer(s) - Name(s) Donna Buckley			
ister	Subject:	510(k) Number K964034/51			-
	То:	The Record - It is my recommendation that the subject	t 510(k) Not	ification:	
		☐ Refused to accept.			
		☐ Requires additional information (other than refuse to	accent)		
	ĺ	☐ Accepted for review	ассори).		
		Is substantially equivalent to marketed devices.			
		☐ NOT substantially equivalent to marketed devices.			
		Other (e.g., exempt by regulation, not a device, dupli	cate, etc.)		
	Is this dev	ice subject to Postmarket Surveillance?	□YES	<b>Z</b> NO	د نی
		ice subject to the Tracking Regulation?	□YES	<b>⊠</b> NO	
		al data necessary to support the review of this 510(k)?	<b>X</b> YES	□NO	
	Is this a pro	escription device?	MYES	□NO	
	Was this 5	10(k) reviewed by a Third Party?	□YES	<b>©</b> NO	
	This 510(k)	) contains:			
	Truthful ar	nd Accurate Statement □Requested ►Enclosed			
	(required	for originals received 3-14-95 and after)			
	□A 510(k)	summary OR A 510(k) statement		i	
	☐ The requ	ired certification and summary for class III devices N/	A	Ÿ	
	The indic	cation for use form (required for originals received 1-1-	96 and after	)	
		er requests under 21 CFR 807.95 (doesn't apply for SE			
	□ No Confi	dentiality     Confidentiality for 90 down to Charles	s):		
		a sample of the second of the			
	Predicate Pro	nduot Code with panel and class: Additional Product	Coéc(s) wit	li panel (optio	ouzi).
-		4 DYB Class II			
j	Review:	Judy Amielson ICOG		Flok-	,
	( <del>Brall</del> Sæ	ich Chief) (Branch Code)		O(Date)	
F	Final Review		no OLT	> Flak	97
	(	(Division Director)	m, The	(Date)	1+
10	levised:11-20-9	76		(vaic)	

# 510(k) "SUBSTANTIAL ROYNVALENCE". Recorder Fold Request #2018-9544, Reteased by CDRH on 09-13-2018 RECORDER OF THE PROPERTY OF THE PROPERTY

er at a GM



N10(%) submissions compare new devices to marketed devices. FDA requests additional intormation if the relationship between marketed and "producesto" (pro-remondments or reclassified post-Amendments) devices is undertail.

- \*\* This decision is normally based on descriptive information alone, but limited testing information is sometimes required.
- \*\*\* Data may be in the \$10(k), other \$10(k)s, the Center's dissilication lifes, or the literature.



## DEPARTMENT OF HEALTH & HUMAN SERVICES

**Public Health Service** Food and Drug Administration

### Memorandum

DATE:

April 11, 1997

TO:

File K964034/S1

FROM:

Donna C. Buckley

Mechanical Engineer, FDA/CDRH/ODE/DCRND/ICDG

SUBJECT:

510(k) Review of Transjugular Access Kit (catheter introducer) -- Angiodynamics

**CONTACT:** Brian Kunst -- (518) 798-1215 ext.190; FAX (518) 798-3625

DECISION: SUBSTANTIALLY EQUIVALENT

74 DYB -- catheter introducer -- 870.1340 -- class II

The sponsor has responded to our January 6, 1997, deficiency letter. In that letter, several questions were asked in order to clarify the device indications and to assess whether the functional characteristics of the device are similar to those of the predicate device.

The device is indicated for percutaneous transjugular access during diagnostic and/or therapeutic interventional procedures. The examples of devices provided in the labeling that are intended for use following liver access are biliary drainage catheters, balloon catheters, and biopsy needles. Since I am unfamiliar with the use of these devices, I requested input from Dr. Daniel Schultz (DRAERD) regarding the device indications and functional characteristics. In response, Dr. Schultz spoke with Anthony Venbrux, Associate Professor of Radiology at Johns Hopkins University, regarding the clinical use of this type of device. Dr. Schultz apprised me of the discussion with Dr. Venbrux and indicated that the stated indications are generally consistent with the current use of the predicate device. For example, an intravascular approach to liver biopsy may be used with patients where a percutaneous liver biopsy would be too hazardous. Also, an intravascular approach may be taken to implant a stent that creates a shunt between the hepatic and portal systems in order to lower the portohepatic pressure gradient. (See attached excerpts and references). Following discussion with several reviewers (and very limited research), I believe that the indications provided in this submission are consistent with the current uses of the predicate.

Since the predicate device, the Cook Rosch-Uchida transjugular liver access set, was declared to be a "modified preamendment" device, no predicate engineering data are on file. Because of the lack of predicate data and the identity of "critical" engineering characteristics of such a device, I believe that the collection of data on the catheter system delivery characteristics (e.g., flexibility, tracking, tensile strength, etc.) and needle characteristics (e.g., sharpness, surface characteristics, mechanical properties, etc.) may not be sufficient to evaluate device safety or efficacy in the absence of both animal studies and clinical data. For this reason, I called Mr. Kunst to discuss my concerns and inquired whether Angiodynamics anticipated collecting animal data or clinical data on the device prior to releasing it to market. Mr. Kunst reported that they had not anticipated collecting animal data in support of the 510(k) submission, however, he agreed to provide me with an article from the Journal of Vascular and Interventional Radiology (fax: 4/24) that reported the clinical use of a device "very similar" to the device under review. The article reported on the use of a particular transjugular liver access set in the placement of Wallstents in 75 patients during transjugular intrahepatic potosystemic shunt (TIPS) procedures. The device referred to in the article was directly compared to the predicate and was reported to "enhance the ability to avoid intervening vital structures and potentially lessen the effect of errant punctures, overcome parenchymal tract resistance in extremely cirrhotic livers, and prevent buckling into the right atrium." (Note: No direct comparisons of clinical data were made to support these conclusions.)

Given the description of the device in the article, I was unable to directly compare the device components to those of the device under review. I called Mr. Kunst and reported that animal data would be essential if he was unable to demonstrate that the devices were either identical or that the differences between them were clinically inconsequential. I received a fax on 4/25 (attached) that compared the components. After review of the supplied information, I believe the device differences to be minor and the referenced clinical information to be supportive of device function. I believe that additional engineering data will not necessarily add new information.

I believe that the prior deficiencies have been adequately addressed and I recommend that this device be found substantially equivalent to 74 DYB (catheter introducer), a class II device.

Down Eng 4/28/9

Donna Buckley

Daniel Schultz, M. D.

(concurrence: evaluation of indications and general function)

Con Ceel Conseelear

the syringe used for anesthetic infiltration, with a small amount of anesthetic left in the syringe. The needle is inserted through a small hole in the skin down to the lesion, which has been visualized. It is then rapidly thrust into the lesion while suction is applied to the syringe, and then it is rapidly withdrawn. Fine-needle aspiration is most often done under ultrasonographic or CT control, so that focal lesions can be sampled. Because the procedure elicits little pain, multiple sticks can be employed, each in a slightly different direction. The material aspirated consists of small fragments of tissue or isolated cells or groups of cells. The aspirate can be smeared on a slide, ejected into appropriate fixatives, or placed in tissue culture media.

Most of the diagnostic interest has been in confirming the presence of a primary hepatocellular carcinoma<sup>51–53</sup> or a hepatoblastoma<sup>54</sup> (Chapter 128). For this purpose, new criteria had to be established for the cytologic diagnosis55,56 and the differentiation of primary carcinoma from metastases in the liver.57 Immunocytochemical methods have also been applied to these specimens to identify the site of origin of a carcinoma; alphafetoprotein and other oncofetal markers, intermediate filaments of the cytoskeleton, and hepatocellular markers are used.58 Special cytologic and cytochemical techniques that can be applied to fine-needle samples are

described later in this chapter. Comparisons have been made between histologic and cytologic methods, and both had about equal diagnostic accuracy; therefore, both are recommended.<sup>59-61</sup> When a carcinoma is suspected but not diagnosed with guided fine-needle biopsy, a core biopsy is recommended to confirm the presence of a benign lesion.62 Hemangiomas can be recognized cytologically<sup>63, 64</sup> and represent a condition that should not indicate a core biopsy. Fine-needle aspirations have also been used to diagnose tumors of the biliary tree and pancreas causing obstructive jaundice. 65, 66 Percutaneous transhepatic catheterization under fluoroscopic control has been suggested as the route to introduce the fine needle to avoid seeding malignant cells along the biopsy track.67 Washings from the needle and from biopsy specimens have been used for cytologic examinations. 68-70 CT guidance<sup>71</sup> has been compared with ultrasonographic guidance, and both seem to be equally effective, although the former was better in the presence of ascites.2

### Transjugular Liver Biopsy 🚁

A completely intravascular approach has been used when a liver biopsy is indicated but the condition of the patient is such that a percutaneous biopsy would be too hazardous.<sup>73-77</sup> The biopsy is best performed by passing a Tru-Cut-style needle inside a 45-cm cardiac catheter into the jugular vein through the superior vena cava directly into a hepatic vein. This must be done under fluoroscopic control and with electrocardiographic monitoring. After the catheter is wedged into a hepatic vein tributary, the needle is thrust into the parenchyma and is triggered via a coaxial cable. Samples

obtained with this Tru-Cut-style needle are almost the same size as percutaneous biopsy specimens. The number of instances in which such a procedure is necessary seems to be small. Often, fine-needle aspiration or laparoscopy is simpler, safer, and less costly in providing the needed information. The femoral vein can also be used for transvenous biopsy with a small biopsy forceps.78

#### Laparoscopic Biopsy

Direct vision is sometimes necessary to sample focal lesions detected by ultrasonography or CT or those that are suspected and not seen by imaging. Laparoscopy, in experienced hands, provides such a technique. It can be used to identify primary<sup>79</sup> or metastatic<sup>80</sup> carcinomas or even the transition of chronic hepatitis to cirrhosis.81 The biopsy is best done percutaneously with a needle, which can be guided by the laparoscopist into the site to be examined. At the same time, a biopsy specimen can be obtained from the uninvolved liver. The Tru-Cut or the Menghini needle can be used, or multiple aspirations with the fine needle can be performed. Morbidity after laparoscopy is minimal. Laparoscopic control has been compared with and preferred to ultrasonographic82 and CT83 control.

#### WHEN NOT TO DO A LIVER BIOPSY

Few contraindications can be listed for liver biopsy. Each experienced operator has a personal set of contraindications, which are often based on personal experience. What follows represents the preferences and prejudices of the authors. The one contraindication most universally agreed on is a bleeding tendency, but the definitions of bleeding tendency vary. A prothrombin time greater than 5 to 6 s from control values is too long for biopsy, although some centers refuse to perform a biopsy if the prothrombin time is 3 to 4 s greater than control values. Postbiopsy bleeding has almost always occurred in patients with normal prothrombin times. Platelet counts less than 50,000 to 60,000/mm<sup>3</sup> are universally considered too low for biopsy, and some centers require more than 75,000 to 80,000/mm<sup>3</sup>. Partial thromboplastin time need not be tested in adults with no history of a bleeding diathesis. The levels of clotting factors can be corrected with fresh frozen plasma and platelet transfusions, but rarely does this have to be done, except for some vascular tumors.84 A biopsy in a patient with greatly disturbed blood coagulation is not likely to yield information not already known and critical to the care of that patient.

Pleural effusion, empyema, pneumothorax, extensive ascites, and peritonitis are reasons to postpone a biopsy until these problems are corrected. Biopsy of cysts and hemangiomas should not be performed blindly. Biopsy of hemangiomas can be performed under direct vision using the laparoscope or with imaging control (see earlier). Extrahepatic biliary obstruction is a relative contraindication, especially if it is long-standing, because



- Frank H, Leodolter I. Praktische Ehrfahrungen mit der ambulanten Leberbiopsie. Wien Klin Wochenschr 1966; 78:756–758.
- Perrault J, McGill DB, Ott DM, et al. Liver biopsy complications in 1000 inpatients and outpatients. Gastroenterology 1978; 74:103–106.
- Sherlock S. Liver biopsy today: the Royal Free Hospital experience. J Hepatol 1984; 1:75–85.
- Maggi G, Alberti A, Bottelli R, et al. La biopsia epatica in regime di day-hospital. Esperienza su una casistica selezionata. Minerva Dietol Gastoenterol 1988; 34:241–244.
- Chawla YK, Ramesh GN, Kaur U, et al. Percutaneous liver biopsy: a safe outpatient procedure. J Gastroenterol Hepatol 1990; 5:94–95.
- Jacobs WH, Goldberg SB, and the Patient Care Committee of the American Gastroenterological Association. Statement on outpatient percutaneous liver biopsy. Dig Dis Sci 1989; 34:322–323.
- tient percutaneous liver biopsy. Dig Dis Sci 1989; 34:322–323.

  34. Brouillette DE, Young-Kul Y, Chien M-C, et al. Use of midazolam for percutaneous liver biopsy. Dig Dis Sci 1989; 34:1553–1558.
- 35. Menghini G. One second biopsy of the liver. Gastroenterology 1958; 35:190–199.
- Rake MO, Murray-Lyon IM, Ansell ID, et al. Improved liver biopsy needle. Lancet 1969; 2:1283.
- 37. Greenwald R, Chirput RO, Schiff ER. Percutaneous aspiration liver biopsy using a large caliber disposable needle: a preliminary report. Am J Dig Dis 1977: 22:1109-1114.
- nary report. Am J Dig Dis 1977; 22:1109–1114.

  38. Mueller PR, Stark DD, Simeone JF, et al. Clinical use of a nonferromagnetic needle for magnetic resonance–guided biopsy. Gastrointest Radiol 1989; 14:61–64.
- Chezmar JL, Keith LL, Nelson RC, et al. Liver transplant biopsies with a biopsy gun. Radiology 1991; 179:447–448.
- Maharaj B, Pillay S. Tru-Cut' needle biopsy of the liver: importance of the correct technique. Postgrad Med J 1991; 67:170–173.
- Letterer H, Drescher T, Busse H-J. Ultraschallgezielte Schneidbiopsie zur Gewinnung von histolgisch auswertbarem Material. Gastroenterol J 1989; 49:156–159.
- Jennings PE, Donald JJ, Coral A, et al. Ultrasound-guided core biopsy. Lancet 1989; 1:1369–1371.
- Colombo M, Del Ninno E, de Freanchis R, et al. Ultrasoundassisted percutaneous liver biopsy: superiority of the Tru-Cut over the Menghini needle for the diagnosis of cirrhosis. Gastroenterology 1988; 95:487–489.
- 44. Solmi L, Muratori R, Brambati M, et al. Comparison between the 21-gauge Urocut needle and the 21-gauge Surecut needle in echo-guided percutaneous biopsy of neoplastic liver lesions. Surg Endosc 1989; 3:38–41.
- Judmaier G, Prior C, Klimpfinger M, et al. Ist der perkutane Leberbiopsie mit der Trucut (Travenol)-Nadel der Menghini-Punktion ueberlegen? Z Gastroenterol 1989; 27:657–661.
- 46. Dani R, de Almeida HR. Estudo prospectivo comparando os resultados obtidos por puncao biopsia hepatica com agulha fina (Chiba) e de Menghini em alguns tumores do figado. Arq Gastroenterol 1987; 24:20–23.
- Farnum JB, Patel PH, Thomas E. The value of Chiba fine-needle aspiration biopsy in the diagnosis of hepatic malignancy: a comparison with Menghini needle biopsy. J Clin Gastroenterol 1989; 11:101–109.
- 48. Parsi B, Guibert JL, Dorcier F, et al. Le diagnostic des lésions hépatiques par ponction à l'aiguille fine. A propos de 407 cas. Sem Hop Paris 1987; 63:2365-2368.
- 49. Buscarini L, Fornari F, Bolondi L, et al. Ultrasound-guided fineneedle biopsy of focal liver lesions: techniques, diagnostic accuracy and complications. A retrospective study on 2091 biopsies. J Hepatol 1990; 11:344-348.
- Fornari F, Civardi G, Cavanna L, et al. Ultrasonically guided fine-needle aspiration biopsy: a highly diagnostic procedure for hepatic tumors. Am J Gastroenterol 1990; 85:1009–1013.
- Trêmolda F, Benevegna L, Drago C, et al. Early detection of hepatocellular carcinoma in patients with cirrhosis by alphafetoprotein, ultrasound and fine-needle biopsy. Hepatogastroenterology 1989; 36:519–521.
- Bru C, Maroto A, Bruix J, et al. Diagnostic accuracy of fineneedle aspiration biopsy in patients with hepatocellular carcinoma. Dig Dis Sci 1989; 34:1763–1769.
- Kung ITM, Chan S-K, Fung K-H. Fine-needle aspiration in hepatocellular carcinoma: combined cytologic and histologic approach. Cancer 1991; 67:673–680.

- 54. Wakely PE, Silverman JF, Geisinger KR, et al. Fine needle biopsy aspiration cytology of hepatoblastoma. Mod Pathol 1990; 3:688–693.
- 55. Sangalli G, Livraghi T, Giorgano F. Fine needle biopsy of hepatocellular carcinoma: improvement in diagnosis by microhistology. Gastroenterology 1989; 96:524–526.
- 56. Cohen MB, Haber MM, Holly EA, et al. Cytologic criteria to distinguish hepatocellular carcinoma from nonneoplastic liver. Am J Clin Pathol 1991; 95:125–130.
- 57. Servoll E, Viste A, Skaarland E, et al. Fine needle aspiration cytology of focal liver lesions: advantages and limitations. Acta Chir Scand 1988; 154:61–63.
- Bedrossian CWM, Davila RM, Merenda G. Immunocytochemical evaluation of liver fine-needle aspirations. Arch Pathol Lab Med 1989; 113:1225–1230.
- Barten M, Kreuzberg M, Kramm G. Vergleich von Feinnadelaspirationszytologie und Histologie beim hepatozelluraeren Karzinom und bei Lebermetastasen. Zentralbl Allg Pathol 1989; 135:147–154.
- 60. Glenthoej A, Sehested M. Histological and cytological fine needle biopsies from focal liver lesions. Intra- and interobserver reproducibility of diagnoses. Acta Pathol Microbiol Immunol Scand C 1989; 97:611–618.
- Glenthoej A, Sehested M, Torp-Pedersen S. Diagnostic reliability of histological and cytological fine needle biopsies from focal liver lesions. Histopathology 1989; 15:375–383.
- Spamer C, Brambs H-J, Kock HK, et al. Benign circumscribed lesions of the liver diagnosed by ultrasonically guided fineneedle biopsy. J Clin Ultrasound 1986; 14:83–88.
- Caturelli É, Rapaccini GL, Sabelli C, et al. Ultrasound-guided fine-needle aspiration biopsy in the diagnosis of hepatic hemangioma. Liver 1986; 6:326–330.
- Nakaizumi A, Iishi H, Yamamoto R, et al. Diagnosis of cavernous hemangioma by fine needle aspiration biopsy under ultrasonic guidance. Gastrointest Radiol 1990; 15:39–42.
- Cohan RH, Illescas FF, Braun SD, et al. Fine needle aspiration biopsy in malignant obstructive jaundice. Gastrointest Radiol 1986: 11:145–150.
- Hall-Crags MA, Lees WR. Fine-needle aspiration biopsy: pancreatic and biliary tumors. AJR 1986; 147:399–403.
- 67. Kuroda C, Yoshioka H, Tokunaga K, et al. Fine-needle aspiration biopsy via percutaneous transhepatic catheterization: technique and clinical results. Gastrointest Radiol 1986; 11:81–84.
- Vogel W, Fagan EA, Bomford A, et al. Wash off liver cytology: a complementary tool to liver biopsy. J Clin Pathol 1986; 39:449–452.
- 69. Axe SR, Erozan YS, Ermatinger SV. Fine-needle aspiration of the liver: a comparison of smear and rinse preparations in the detection of cancer. Am J Clin Pathol 1986; 86:281–285.
- Fagelman D, Chess Q. Nonaspiration fine-needle cytology of the liver: a new technique for obtaining diagnostic samples. AJR 1990; 155:1217–1219.
- 71. Lüning M, Schroeder K, Wolff H, et al. Percutaneous biopsy of the liver. Cardiovasc Intervent Radiol 1991; 14:40–42.
- Murphy FB, Barefield KP, Steinberg HV, et al. CT- or sonography-guided biopsy of the liver in the presence of ascites: frequency of complications. AJR 1988; 151:485–486.
- Roesch J, Lakin PC, Antonovic R, et al. Transjugular approach to liver biopsy and transhepatic cholangiography. N Engl J Med 1973; 289:227–231.
- 74. Gamble P, Colapinto RF, Stronell RD, et al. Transjugular liver biopsy: a review of 461 biopsies. Radiology 1985; 157:589–593.
- 75. Lebrec D. Transvenous (transjugular) liver biopsy. Intern Med Specialist 1987; 8:81–94.
- Steadman C, Nathan N, Teague C, et al. Transjugular liver biopsy—an Australian experience. Aust N Z J Med 1988; 18:836– 840.
- 77. Vlavianos P, Bird G, Portmann B, et al. Transjugular liver biopsy: use in a selected high risk population. Eur J Gastroenterol Hepatol 1991; 3:469–472.
- 78. Mewissen MW, Lipchik EO, Schreiber ER, et al. Liver biopsy through the femoral vein. Radiology 1988; 169:842–843.
- Jeffers L, Spieglman G, Reddy R, et al. Laparoscopically directed fine needle aspiration for the diagnosis of hepatocellular carcinoma: a safe and accurate technique. Gastrointest Endosc 1988; 34:235–237.

The combination of sclerotherapy with propranolol has been used to reduce the rebleeding rate related to longterm injection sclerotherapy.391 Some reported a decrease in the frequency of rebleeding without an increase in complications before obliteration.392, 393 This combination appeared to be most beneficial in patients with advanced alcoholic liver disease. Others found an increased complication rate (encephalopathy, difficult resuscitation) without reduction of the frequency of rebleeding or its severity.391

#### Surgical Therapy

When patients have bled from varices, surgery is a third therapeutic option. The history of portal decompressive surgery spans 4 decades of the 20th century. In the early 1990s it is not considered a first line of treatment, but in the 1970s therapeutic portacaval anastomosis was the mainstay of therapy. Four randomized trials examined the efficacy of central shunts (end-toside portacaval anastomosis was done in more than 90% of the patients) in preventing recurrent hemorrhage in patients with cirrhosis (the majority of alcoholic etiology). Control groups received no treatment.

The benefits of shunt surgery over conventional treatment were marginal in a total of 292 patients.340,341, 394, 395 Surgical mortality ranged from 7% to 21%, and 5% of patients had a compromised shunt. Survival was slightly increased in the shunt group;340, 341, 395 it was shortened in a fourth study in 3 of the trials.394 As expected, the incidence of bleeding was reduced in the surgical group but encephalopathy was increased. When the studies are combined, the survival curve was

still slightly better in the surgical group.

The negative experience with postshunt encephalopathy fueled interest in developing new selective decompressive procedures.329 The distal splenorenal shunt is technically more demanding, requiring an arduous dissection. Complications arise, in part, from an initial learning phase. Seven trials compared this selective shunt to central shunts in the prevention of recurrent hemorrhage. Table 102-6 lists these studies.396-402 The majority of the patients had alcoholic cirrhosis (78%). Surgical mortality was higher in the distal splenorenal group, reflecting the previously mentioned difficulties; long-term survival was not drastically different. The frequency of postshunt variceal hemorrhage was disturbingly high in some studies, with a high prevalence of shunt occlusion. The prevalence of postshunt encephalopathy, however, was markedly different among the trials. Whereas 3 studies showed a lower rate with distal splenorenal shunting, 3 others did not. The tools used to define encephalopathy varied among studies, but in one report that examined quality of life and hospitalizations for encephalopathy, 398 a clear benefit for the selective procedure was noted.

Patients with nonalcoholic cirrhosis may have a better response to the distal splenorenal shunt,403 although this was not confirmed by others. 402 The progression of liver disease may result in the development of collaterals from the mesenteric bed and in a loss of selectivity,

the original rationale for the operation.<sup>404</sup> Claims for better preservation of liver function after distal splenorenal shunts405 require further confirmation.

Modifications of the procedure include a pancreaticosplenic venous disconnection406 to avoid a new pathway for collateralization (Fig. 102-9). The dissection, however, can be particularly difficult in alcoholic patients with chronic pancreatitis. Distal splenorenal shunts are not suitable for patients with tense ascites; transient accumulation of ascitic fluid during volume resuscitation, however, is common and should not preclude the procedure if hepatopetal flow is demonstrated preoperatively. The procedure must be performed by an experienced surgeon; the availability of alternative methods of treatment has reduced the number of shunt procedures worldwide.

Postshunt encephalopathy continues to be a major concern when recommending surgical decompression. No methods predict the development of postsurgical changes in mental state. Transient encephalopathy that occurs around the bleeding episode does not preclude a shunt procedure.407 Preoperative liver function does not predict the neurologic response to the procedure. 408 Experience with the small-diameter portacaval H graft,409 in which a preserved portacaval pressure gradient and higher hepatopetal flow result in lower rates of encephalopathy, suggests that a key factor in the incidence of postshunt alteration in mental state is the rate of delivery of neurotoxins from the splanchnic bed to the systemic circulation.410

Three studies compared distal splenorenal shunts with sclerotherapy in alcoholic cirrhotic patients (80% to 90% of the study population). $^{411-413}$  The rebleeding rate was higher after sclerotherapy, but no differences in survival were seen in all studies. Of patients in the sclerotherapy group, 6% to 31% ultimately underwent shunt surgery. Hepatic portal perfusion was maintained better by sclerotherapy than by selective shunting. One study, comparing end-to-side portacaval anastomosis to sclerotherapy, found surprisingly similar incidences of encephalopathy.367

Surgery is proposed for patients considered unresponsive to sclerotherapy. However, sclerotherapy requires a long-term commitment from the patient. The choice of therapy should be based on the assessment of patient's compliance, geographic considerations, and

the availability of local expertise.

### **Future Directions**

### Transjugular Intrahepatic Portosystemic Shunt Stent

A radiologic technique for the percutaneous placement of an intrahepatic shunt has been developed.415 The procedure is rapid, requiring 2 to 3 hours to complete, and is performed without anesthesia. It involves localization of the portal vein by ultrasonography and transjugular establishment of a parenchymal tract between a central hepatic and portal vein (Fig. 102-10). An expandable wire mesh stent (diameter 8 to 12 mm)

Table 102-6. RANDOMIZED CONTROLLED TRIALS COMPARING DISTAL SPLENORENAL SHUNTS VERSUS NONSELECTIVE SHUNTS IN PREVENTING RECURRENT VARICEAL HEMORRHAGE\*

		Number of Patients			Operative Mortality (%)		Variceal Hemorrhage		Shunt Occlusion (%)		PSE (%)	
Trial	Control Shunt	PSS	DSRS	Cirrhosis (%)	PSS	DSRS	PSS	DSRS	PSS	DSRS	PSS	DSRS
Millikan et al. <sup>396</sup>	RCS,	29	26	72	10	12	31	22	29	14	76	26
Reichle et al. <sup>397</sup> Langer et al. <sup>398</sup> Grace et al. <sup>399</sup> Harley et al. <sup>400</sup> Fischer and McKinley <sup>401</sup> Spina et al. <sup>402</sup> Totals	MCS MCS PCA PCA PCA	14 40 38 27 27 47 47	13 38 43 26 29 46 221	93 82 90 100 84 46 78	7 0 13 7 7 2 7	8 13 9 12 10 2 9	5 12 4 7 2 10	1R† 9 18 30 11 9 16	15 N 12	IR IR 17 IR 23 IR	48 45 32 16 17 39	36 51 39 8 21 30


<sup>\*</sup>PSS = portosystemic shunt; DSRS = distal splenorenal shunt; PSE = portosystemic encephalopathy; RCS = renocaval shunt; MCS = mesocaval shunt; PCA = portacaval anastomosis; PSRS = proximal splenorenal shunt.

tNR = not reported.

Adapted from Grace ND, Conn HO, Resnick RH, et al. Distal splenorenal vs portal-systemic shunts after hemorrhage from varices: a randomized controlled trial. Hepatology 1988; 8:1475–1481.

is placed and a balloon expanded to achieve a reduction of the portal venous pressure gradient to less than 12 mm Hg. A neointimal lining of the stent forms within a few weeks in the dog model.<sup>415</sup>

Experience in nonsurgical candidates with portal hypertension and refractory variceal bleeding<sup>416</sup>. <sup>417</sup> showed that TIPS lowered the mean portohepatic pressure gradient with cessation of bleeding. When patients required a greater portosystemic pressure gradient reduction, the stent was further expanded by percutaneous balloon angioplasty. <sup>417</sup> Followup endoscopies within 1 week showed a reduction in size of the varices in all patients. All shunts were patent within the mean followup duration of 5 months, although additional balloon angioplasty procedures were needed for minor strictures in 2 of 8 patients.



**Figure 102–9.** Development of new collaterals after distal splenorenal shunt results in loss of selectivity of the operation; transgastric and pancreatic collaterals shunt blood from the portal vein to the splenic venous drainage. (From Warren WD, Millikan WJ, Henderson JM, et al. Splenopancreatic disconnection. Improved selectivity of distal splenorenal shunt. Ann Surg 1986; 204:346–355.)

Further studies are needed to compare present therapeutic modalities with TIPS in terms of rebleeding rate, complication rates, and long-term patency. TIPS was highly effective in controlling both the acute and recurrent hemorrhages in a small uncontrolled prospective study of cirrhotic patients with advanced Child's class and refractory variceal bleeding. 135 Portal pressures were reduced from a mean of 34 to 22 mm Hg, but portal hepatic gradients were not measured. No procedural or intrahepatic complications were reported and hepatic function was preserved. Potential complications, however, include stent migration, hemobilia, and arteriovenous fistula. A major advantage of this procedure is the preservation of the porta hepatis for future transplantation. Results on long-term survival are awaited.

### Liver Transplantation

Long-term survival rates with liver transplantation have been reported from some centers to be as high as 70% at 3 years. This approaches the survival seen after elective shunt surgery and allows consideration of recurrent variceal bleeding as an indication for transplantation in patients with decompensated liver disease. However, organ supply and health care costs prohibit its widespread use, and transplantation should be reserved for individuals with variceal hemorrhage and concomitant hepatocellular failure. Still, hemorrhage should be controlled by other means before considering transplantation.

### ADDITIONAL BLEEDING SITES

### Congestive Gastropathy

Endoscopic descriptions of portal hypertensive gastropathy include a mosaic pattern of 2 to 6 mm erythematous patches separated by a fine yellow or white lattice or, in a more severe form, a granular mucosa with cherry red spots, also termed scarlatiniform pat-

Records processed under FOIA Request #2018-3544; Released by CDRH on 09-13-2018

# ANGIODYNAMICS, Inc.

Bill To:

603 Queensbury Avenue Queensbury, NY 12804 Ship To:

543 Queensbury Avenue Suite 4

Queensbury, NY 12804

Phone Fax 518-798-1215 518-798-3625

FAX TRANSMITTAL

TO:

EDA-ODR

ATTN:

r-over M

CC:

FROM:

JET MANAJion

SUBJECT:

1

(Including this cover sheet)

Message:

PONNA

TE you have any other questions regarding this submission or the

Brown Killing Line Det A CALL.

TOWAY & TOMORION (5/9-5/6)

THE INFORMATION CONTAINED IN THIS TRANSMISSION IS PRIVILEGED AND CONFIDENTIAL .IT IS INTENDED ONLY FOR THE USE OF THE INDIVIDUAL OR ENTITY NAMED ABOVE. IF THE READER OF THIS MESSAGE IS NOT THE INTENDED RECIPIENT, YOU ARE REFERY NOTIFIED THAT ANY DISSEMINATION, DISTRIBUTION OR COPYING OF THIS COMMUNICATION IS STRICTLY PROHIBITED. IF YOU HAVE RECEIVED THIS COMMUNICATION IN ERROR, PLEASE NOTIFY US IMMEDIATELY BY TELEPHONE AND RETURN THE ORIGINAL MESSAGE TO US AT THE ABOVE ADDRESS VIA THE U.S. POSTAL SERVICE, THANK YOU.

Records processed under FOIA Request #2018-3544; Released by CDRH on 09-13-2018

Pa	ige	1	of	1

510(k) Number (if known): K964034

Device Name: Transjugular Access Set

Indications For Use:

(Per 21 CFR 801.109)

The AngioDynamics, Inc. Transjugular Liver Access Set is intended to be used for transjugular liver access during diagnostic and interventional procedures.

NEEDED)			
		s in independent of the second	
Concurrence of	CDRH, Office of L	Device Evaluation (ODE)	
Prescription Use	or	Over-The-Counter Use	

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF

(Optional Format 1-2-96)

Records processed under FOIA Request #2018-3544; Released by CDRH on 09-13-2018

# AngioDynamics<sup>®</sup>

INCORPORATED

April 25, 1997

Donna Buckley
U.S. Food and Drug Administration
Center for Devices and Radiological Health
9200 Corporate Drive
Rockville, Maryland 20850

Re:

Additional Information Request

K964034 - ANGIOLYNAMICS® Transjugular Access Set

Dear Ms. Buckley:

The following information is provided in response to your request for additional information regarding the above referenced premarket notification.

The 21 gauge needle is manufactured by Needletech, North Attelboro, Massachussetts. The 14 gauge needle is manufactured by Hawk Medical, Gainesville, Florida.

The following table compares the Access Set described in K964034 to the device described by Hawkins in the article published in the Journal of Vascular and Interventional Radiology.

Device presented in K964034	Device presented in JVIR article Nov. 1994
8 French introducer sheath & hemostasis valve	12 French guiding catheter & hemostasis valve
8 French catheter	8 French catheter
5 French AngiOptic Catheter	5 French catheter
14 gauge curved needle with "extra blunt" tip	14 gauge curved needle with blunt tip
21 gauge beveled tip needle	22 gauge beveled tip needle

#### 8 French introducer sheath vs. 12 French guiding catheter:

The 8 French introducer sheath presented in K964034 and the 12 French guiding catheter presented in the JVIR article are the same device. The 8F refers to the ID or the maximum diameter of the catheter that can be introduced through the sheath. The OD is 12F. Since industry standard for introducers is to label them with the ID or the maximum size of the device it is intended to introduce, it will be referred to as an 8F introducer sheath.

#### 8 French PTFE catheter vs. 8 French catheter:

The 8 French PTFE catheter presented in K964034 and the 8 French catheter presented in the JVIR article are the same device.

603 Oneensbury Ave. - Queensbury, NY 12804 USA • 518-798-1215 \* 800-77ANGIO \* fax: 518-798-3625 E-mail: http://www.angiodynamics.com



### 5 French AngiOptic catheter vs. 5 French catheter.

1-1F-17 11 12 EEV

The 5 French catheter used in device described in the JVIR article was constructed from a Teflon shaft, while the device presented in K964034 uses a commercially available nylon shaft straight AngiOptic catheter (K954423). The reason for this change was that the straight AngiOptic is already in our product line and the purpose of the 5F catheter is for contrast injection — consistent with the design goal of the AngiOptic. There are no dimensional differences between the two devices.

# 14 gauge curved needle with "extra-blunt" tip vs. 14 gauge curved needle with blunt tip:

As figure 1 indicates the change from the blunt tip to the "extra-blunt" tip was to reduce the possibility of severing the tip of the 5 French catheter while retracting it through the 14 gauge curved needle.




Figure 1 - Comparison of two styles of curved needle

### 21 gauge needle vs. 22 gauge needle:

The reason for increasing the size of the beveled point needle from 22 gauge to 21 gauge was to allow the use of an .018" guidewire instead of an .016" guidewire. A .018" can in some instances be easier to track through the vasculature and gives the physician more flexibility.

The differences between the device proposed in the 510(k) and the device described in the article are very slight and would not affect the product's safe and effective use.

If there are any further questions regarding this application, please feel free to contact me at 518-798-1215 extension 123.

Sincerely,

ANGIODYNAMICS, Inc.

Brian Kunst

Director, Regulatory Affairs and Quality Assurance



		YNAMICS, Inc.				
	Phono: 518-798-121	nue, Queensbury, New York 12804 15 Fax: 518-798-3625				
FAX		Date: April 21, 1997  Number of pages including cover sheet: 5				
To:	Donna Buckley Division of Cardiovascular Device	From: Brian S. Kunst Director, Regulatory Affairs and Quality Assurance				
Phone: Fax phone: CC;	301-480-4204	Phone: 518-798-1215, ext. 123  Fax: 518-798-3625				
REMARKS: Urgent For your review Reply ASAP Please comment  Re: K964034 Additional Information						
Dear Ms. Buckley,						
Hawkins specified the experience as descri	be features of the AngioDynamics Transjugula	Portal Venous Access System by Irvin Hawkins, MD. Dr. r Access Set described in K964034. based on his clinical nee that the AngioDynamics Transjugular Access Set will be				
If you have any que:	tions, please call me at 518-798-1215, extensi	on -33				
Sincerely,		<b>]</b>				
ANGIODYNAMIC	S <sup>®</sup> , Inc.					
Brian Kunsı Director of Regulate	ry Affairs/Quality Assurance					

# Technical Developments and Instrumentation

# Fine-Needle Transjugular Portal Venous Access System<sup>1</sup>

Scott R. Kerns, MD Frank W. Sabatelli, MD<sup>2</sup> Irvin F. Hawkins, Jr, MD

THE transjugular intrahepatic portosystemic shunt (TIPS) has become a viable option in the treatment of patients with portal hypertension. Numerous series have shown it effectively lowers portal pressures, thereby reducing the chance of variceal hemorrhage and allowing resorption of assites (1-4). Those results are achieved with a lower procedural morbidity and mortality compared with those associated with surgical shunting. Excluding hepatic encephaiopathy and late shunt stenosis or occlusion, the complication rate of TIPS procedures has been reported to be less than 10% (5). Although rare, serious and potentially fatal complications may occur including hepatic artery occlusion or pseudoaneurysm formation, arterioportal or arteriovenous fistula formation, hemobilia, and intraperitoneal hemorrhage. Various techniques and equipment for transjugular portal venous access have been proposed in part to decrease procedural complications while maintaining a high technical success rate (1,6,7). With these goals in mind we have developed and clinically tested a modified transjugular portal vonous access system, and in this report describe our experience using it to perform TIPS procedures in 75 patients

PATIENTS AND METHODS

The fine-needle set is a coaxial system consisting primarily of five components (AngioDynamics, Glens Falls, NY) (Figs 1, 2). The inner component consists of an 82-cm-long, 22-gauge, thin-wall puncture needle inside a 70-cm-long, 5-F catheter. A 0.016-inch nitinol guide wire with a highly radiopaque tip (Microvena. White Bear Lake, Minn) is placed into the needle. When the needle is withdrawn into the 5-F catheter, a stop on the shaft of the needle exactly aligns their tips. This prevents accidental advancement of the needle through the wall of the catheter.

The outer component consists of a

blunt-tip, 60-cm-long, 14-gauge metal cannula inside a 57-cm-long, 6-F catheter. The cannula has a 30° curve over its distal 4 cm that directs the course of the emaller puncture needle. A metal handle on the hub of the cannula indicates the direction of its curved tip and acts as a handle for torquing the cannula. The entire assembly is placed through a 44-cm-long, 12-F, thick-walled Teflon guiding catheter with a hemostasis valve and side port.

This fine-needle system was used to perform TIPS in 75 patients between May 1991 and March 1994. The most common presentation was recurrent variceal hemorrhage refractory to sclerotherapy. After the procedure, patients were observed at least overnight with daily scrial measurements of hematocrit levels. A color Doppler flow imaging examination was performed the morning after the procedure to confirm shunt patency and to evaluate for vessel occlusion or pseudoaneurysm, tiliary obstruction, or any change in the amount or character of ascites.

The general TIPS technique has been previously described (1-3); however, several modifications were used. The right internal jugular vein was entered under direct ultrasound (US) guidance with a 21-gauge needle. Once in the hopatic vein, the 12-F guiding catheter was retracted and the curved metal cannula rotated toward the expected location of the right portal vein. The 22-gauge needle was advanced in 1-2-mm increments, and the area was tested by gently probing with the preloaded 0.016 inch guide wire or injecting contrast material via a Thuby-Borst adapter. This was continued until the portal vein was entered, un orrant puncture of another structure occurred, or the needle had been advanced 5 or 6 cm. Errant punctures included puncture of the hepatic artery, bile duct. or liver capsule. If an errant puncture occurred or the needle advanced 5 or 6 cm without entering the purtal vein, the needle was retracted, the cannula redirected, and a new sories of passes per-

Index terms: Catheters and catheterization, technology, 957.123 Portal vein, 957 123 Shunts, portosystomi 959.463

JVIR 1994; 5:835-8:7

1 From the Department of Radiology, Shands Teaching Hospital and Clinies, University of Florida, 1600 SW Archer Rd. Gainesville, FL 92010. From the 1994 SCVIR annual meeting. Received March 22, 1994 revision requested June 1; revision received June 20; accepted June 24. Address reprint requests to S.R.K.


2 Corrept address: CPO Radiology Associates.
Plane, Tox.


© SCVIR, 1994

835

## 836 • Journal of Vascular and Interventional Radiology

November December 1994





Figures 1, 2. (1) Desassembled fine-needle transjugular portal venous access system. From top to bottom: 22-guage needle with stop on proximal shaft, 5-F catheter, 14-gauge blunt metal cannula, 8-F catheter, and 12-F guiding catheter with homostasis valve and side port. (2) Assembled fine-needle portal venous access system.

Once the 22-gauge needle had entered the portal venous system without traversing the hepatic artery or biliary system, the 0.016-inch guide wire was advanced followed by the needle and 5-F catheter. Portography was performed to confirm that the right portal vein had been entered. The 8-F and 12-F catheters were then advanced into the portal vein over a guide wire followed by balloon dilation and placement of a Wallsten (Schneider, Minneapolis, Minn). The sturdy 12-F catheter provided direct access to the main portal vein for balloon expansion and stent placement without buckling or resistance.

#### RESULTS

We were able to place a functioning shunt in 72 of 75 patients, for a technical success rate of 96%. Among our failures, one patient had previous portal vein thrombosis with partial recanalization on imaging studies. In this patient, we were unable to enter a recanalized branch despite approximately 50 passes. In another patient the posterior branch of the right portal vein thrombosed along with the abunt after dilation of the shunt. The patient returned later for

successful placement of a shunt to the left portal vein. The third failure occurred in a patient with a large liver and very small portal vein branches. We were unable to enter the portal venous system, and he later underwont successful creation of a distal splenorenal shunt.

2

We had no difficulty advancing the 0.016-inch guide wire into the portal system, and it never was sheared off by the sharp needle. However, we did not attempt to pass the needle around acute turns for fear of this occurring.

No patients developed clinical signs of hemorrhage or a drop in hematocrit, and no evidence of vascular injury, biliary obstruction, or hemoperitoneum was seen on the postprocedure US examinations. The hepatic artery, biliary system, and liver capsula were inadvertently punctured in approximately 5%, 10%, and 20% of patients, respectively.

No extant puncture of the gallbladder, colon, or right kidney occurred. In two patients (3%) the main portal vein, which is believed to be extrahepatic, was entered without sequelae. In all of these patients new passes were performed to complete the TIPS without violating adjacent structures.

We found that in patients with small livers, the portal bifurcation was frequently almost directly anterior to the hepatic vein, which necessitated a sharp turn anterior to make a successful puncture. In these cases, the angle of the blunt cannula can be increased by carefully bending it by hand to the desired amount. Also, the 22-gauge needle may be bent as well, although this should be done with the 0.016 inch guide wire in side it to avoid compromising the needle lumen. The curve of the hent 22-gauge needle naturally follows and extends the curve of the metal cannula.

#### DISCUSSION

In general, TIPS is considered a relatively safe means of portal decompression in patients with portal hypertension. Although most procedural complications either are treated conservatively or have no significant impact on the patient's care, serious and potentially fatal complications can and do occur (5,5-10). To improve the safety of this procedure, thereby potentially increasing its clinical usefulness, we have attempted to refine what may be the most difficult and dangerous part of this procedure the transjugular portal venous

The fine-needle access system described is a modification of other currently available systems used for TIPS (1,11,12). Unlike the system described Rösch et al (12), we used a small-gauge hollow needle instead of a trocar to make passes toward the portal vein. Passes are short, and the location of the needle tip is checked before proceeding further. This enables one to avoid traversing or transecting other vital structures in the vicinity of the portal vein. Most investigators report using a larger gauge sharp Colapinto needle or the Rüsch set with its trocar stylet to perform the puncture. Thrusts are made to-ward the portal vein, and aspiration is applied while the needle or catheter is withdrawn until blood return is obtained. Potentially, a vital structure could be entered on the way to the por tal vein without the operator knowing it. In approximately one-half of patients another vital structure lies in the path between the right|hepatic vein and the portal bifurcation (13). This may result in arterial injury, fistula formation, or communication of the shunt with the biliary tree. In addition to procedural morbidity, the latter should probable be

avoided since preliminary data show this may incite pseudointimal hyperplasia leading to shunt failure (14). With this system the tract is tested before the sheath is advanced and the shunt is placed, thereby avoiding these problems.

The guiding catheter, which is longer and has a thicker wall than the thinwalled sheath provided with the Rosch set, provides access to the portal vein in all patients without buckling into the right atrium. We and others have found that with the Rosch set this has not always been the case (15). Its length allows it to be placed into the portal vein completely eliminating any resistance in the parenchymal tract to the balloon catheter or Wallstent deployment device


In conclusion, we believe this fineneedle transjugular portal venous access system can be used safely and effectively for the TIPS placement. This system enhances the ability to avoid intervening vital structures and potentially lessen the effect of errant punctures, overcome parenchymal tract resistance in extremely cirrhotic livers, and prevent buckling into the right atrium. Ultimately, the choice of access system will dopend on how much significance the operator places on these factors.

#### References


- Ring EJ, Lake JR, Roborts JP, et al Using transjugnlar intrahepatic porto systemic shunts to control arriceal bleeding before liver transplantation. Ann Intern Med 1992; 116:304-309.
   LaBerge JM, Ring EJ, Cordon RL, et al.
- LaBerge JM, Ring EJ, Cordon RL, et al. Creation of transjugular intrahepatic portosystemic alimits with the Wallstent endoprosthesis: results in 100 patients. Radiology 1993; 187:413-420.
- Rousseen H, Vinel J, Bilbao JI, et al.
  Transjugular intrahepatic portosystemic
  shunts using the Wallstent prosthesis: a
  follow-up study. Cardiovasc Intervent
  Radiol 1994; 17:7-11
- 4 Ferral H. Bjarnason H. Wegryn SA, et al. Refractory ascites and experience

- in treatment with transjugular intrahepatic purlosystemic shunt. Radiology 1993; 189:790-801.
- Freedman AM, Sanyal AJ, Tisnado J, et al Complications of transjugular intrahepatic portosystemic shunt: a compreheusive review. RadioGraphico 1993; 13:1185-1210.
- Maynar M, Cabrera J, Pulido-Duque JM, et al. Transjugular intrahepatic porto systemic shunt: early experience with a flexible trocar/catheter system. AJR 1993; 161:301-306.
- 7 Rozenbut G, Dei Guercio LRIV. Combined transmesenteric and transjugular approach for intrahepatic portosystemic shupt placement. JVIR 1993; 4:661-666.
- abunt placement. JVIR 1993; 4:661-666.

  McGse AL, Reas CR, Niblett RL, Diamond NG, Lee SP. Complications after TIPS: predictive variables (abstr). JVIK 1994; 5:3.
- Rocele M. Puncture of the portal bifurcation, a fatal complication of TIPS (let ter), RadioGraphics 1993; 13:1184.
- 10 Haskal ZJ, Pentecost NJ, Rubin RA. Hapatic arterial injury after transjugular intrahepstic portosystemic shunt placement; report of two cases. Radiology 1993; 188:85-88.
- Uchida BT, Putnam JS, Rosch J. "Arraumatic" transjugular needle for portal vein puncture in swine. Radiology 1987; 183:580-581.
- Rosch J, Uchida DT, Darton RE, Keller FS. Coaxial catheter-needle system for transjugular portal vein entrance. JVIR 1993; 4:145-147.
- 13 Uflocker R, Roichort P, D'Albuquerque LC, Silva AO. Liva anatomy applied to the placement of transjugular intrahepatic portosystemic shunts. Radiology 1994; 191:705-712.
- LuBerye JM, Ferrell LD, Ring EJ, Cordon RL. Histopathologic study of sto notic and occluded transjugular intrahepatic portosystemic shunts. JVIR 1993; 4:779-786.
- McKusick MA, Johnson CM, Williams IIS. Shiffening connuls to facilitate tract dilation for transjunds: introhepatic portosystemic shunt placement. JVIR 1993: 4:787-788.



Previous fine needle system. Thin-walled, "blunt" 14gauge Ross needle "slices" 5 French catheter during retraction. .-[1-]\* (C:20.1 ::0 :



New "extra blunt" Ross needle has a rounded tip which will not cut inner 5 French catheter. The coaxial 8 French catheter has been eliminated. Note: The 21-gauge needle is advanced during search for portal vein. The 5 French catheter traverses liver chly after 21-gauge needle enters contail vein.

Food and Drug Administration Center for Devices and Radiological Health Office of Device Evaluation Document Mail Center (HFZ-401) 9200 Corporate Blvd. Rockville, Maryland 20850

February 05, 1997

ANGIODYNAMICS, INC. 603 QUEENSBURY AVE. QUEENSBURY, NY 12804 ATTN: BRIAN KUNST 510(k) Number: K964034
Product: TRANSJUGULAR
ACCESS SET

The additional information you have submitted has been received.

We will notify you when the processing of this submission has been completed or if any additional information is required. Please remember that all correspondence concerning your submission MUST be sent to the Document Mail Center (HFZ-401) at the above letterhead address. Correspondence sent to any address other than the one above will not be considered as part of your official premarket notification submission. Because of equipment and personnel limitations we cannot accept telefaxed material as part of your official premarket notification submission, unless specifically requested of you by an FDA official.

The Safe Medical Devices Act of 1990, signed on November 28, states that you may not place this device into commercial distribution until you receive a letter from FDA allowing you to do so. As in the past, we intend to complete our review as quickly as possible. Generally we do so 90 days. However, the complexity of a submission or a requirement for additional information may occasionally cause the review to extend beyond 90 days. Thus, if you have not received a written decision or been contacted within 90 days of our receipt date you may want to check with FDA to determine the status of your submission.

If you have procedural or policy questions, please contact the Division of Small Manufacturers Assistance at (301) 443-6597 or at their toll-free number (800) 638-2041, or contact me at (301) 594-1190.

Sincerely yours,

Marjorie Shulman Supervisory Consumer Safety Officer Premarket Notification Section Office of Device Evaluation Center for Devices and Radiological Health Records processed under FOIA Request #2018-3544; Released by CDRH on 09-13-2018

February 4, 1997

U.S. Food and Drug Administration Center for Devices and Radiological Health Document Mail Center (HFZ-401) 9200 Corporate Drive Rockville, Maryland 20850



Re: Additional Information request K964034 - AngioDynamics® Transjugular Access Set

Dear Document Control Staff:

The following information is provided in response to your request for additional information, dated January 6, 1997, regarding the above referenced premarket notification.

### 1). Indications For Use

ANGIODYNAMICS, Inc. wishes to revise the indication for use for the Transjugular Access Set. The revised indication for use is as follows:

"Transjugular liver access in diagnostic and interventional procedures"

This is identical to the labeled indication for use of the predicate device. Attachment 1 contains the revised instructions for use.

### 2). Beveled Needle

The purpose for using a hollow needle in the AngioDynamics® Transjugular Access Set instead of the solid pencil point needle used in the Cook® Rosch-Uchida Transjugular Liver Access Set is to enable contrast injections through the hollow needle during the procedure to assure proper device placement. The Transjugular Access Set was developed in cooperation with Dr. Irvin Hawkins, Professor of Radiology and Chief of Interventional Radiology at Shands Teaching Hospital, Gainesville, FL, whose clinical experience guided the design of this feature. Note that the O.D. of the pencil point needle is 20ga (.035"), which is equivalent to the AngioDynamics, Inc. needle 21ga (.032").

Questions? Contact FDA/CDRH/OCE/DID at CDRH-FOISTATUS@fda.hhs.gov or 301-796-8118

Ø

Both the AngioDynamics, Inc. Transjugular Access Set and the Cook® Rosch-Uchida Transjugular Liver Access Set contain a curved component used for directing the catheter. Both are 14ga in diameter and are made of stainless steel. AngioDynamics, Inc. refers to the shape of the needle as a Ross needle, while Cook® refers to theirs as a stiffening cannula. There are no functional differences.

### 3). Package Label

Originally, the needle that was intended to be used in the Transjugular Access Set was a 22g thin walled needle (.0280 O.D. x .0205 I.D.). A slight variation was made to this by changing the 22g thin walled needle to a 21g thin walled needle (.032" O.D. x .024" I.D.). The purpose of this change was to enable the use of an .018" guidewire, since the 22g needle limited guidewire size to .016". A revised package label is included as Attachment 2.

If there are any further questions regarding this application, please feel free to contact me at 518-798-1215 extension 190.


Sincerely, AngioDynamics, Inc.


Jeffrey R. Mannion Regulatory Affairs Associate





# Attachment 1












# Attachment 2





Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

## JAN 6 1997

Mr. Jeffrey R. Mannion Regulatory Affairs Associate Angiodynamics, Incorporated 603 Queensbury Avenue Queensbury, New York 12804

Re: K964034

Transjugular Access Set Dated: September 30, 1996 Received: October 8, 1996

Dear Mr. Mannion:

We have reviewed your Section 510(k) notification of intent to market the device referenced above. We cannot determine if the device is substantially equivalent to a legally marketed predicate device based solely on the information you provided. To complete the review of your submission, we require the following additional information:

- Please address the following issues regarding the indications for use:
  - a. Each indication for use needs to be specified and supported. You indicate in Section 5 that your device is "intended to be used for transjugular access in diagnostic and interventional procedures, such as biliary drainage, liver access and biopsy, placement of vena cava filters, and cholangiography." However, the device to which you claim substantial equivalence is indicated only for liver access. Therefore, please address the following: (1) revise the indication for your device to specify the areas where it is to be used; (2) identify a predicate device for each indication; and (3) demonstrate that your device is substantially equivalent to the predicate by providing comparative data and performance data.
  - b. Please list the vena cava filters that you intend to place with your device and provide either an explanation or testing to demonstrate the adequacy of your device to place these filters without incident.



### Page 2 - Mr. Jeffrey R. Mannion

- 2. The predicate device includes a pencil point needle, whereas your device includes a Ross needle and a beveled tip needle. Please provide either the trade names and manufacturers of devices that use the 14 gauge curved Ross needle and the 21 gauge beveled tip needle for the same indications or provide an article from a peer reviewed journal that indicates that these needle tips and gauges are adequate for the stated indications.
- 3. The package label indicates that the kit includes a 22 gauge beveled point needle, whereas the rest of the file references a 21 gauge beveled point needle. Please clarify this discrepancy.

We believe that this information is necessary for us to determine whether or not this device is substantially equivalent to a legally marketed predicate device with regard to its safety and effectiveness.

You may not market this device until you have provided adequate information described above and required by 21 CFR 807.87(k), and you have received a letter from FDA allowing you to do so. If you market the device without conforming to these requirements, you will be in violation of the Federal Food, Drug, and Cosmetic Act (Act). You may, however, distribute this device for investigational purposes to obtain clinical data if needed to establish substantial equivalence. Clinical investigations of this device must be conducted in accordance with the investigational device exemption (IDE) regulations.

If the information, or a request for an extension of time, is not received within 30 days, we will consider your premarket notification to be withdrawn and your submission will be deleted from our system. If you submit the requested information after 30 days it will be considered and processed as a new 510(k); therefore, all information previously submitted must be resubmitted so that your new 510(k) is complete.

Bil

Page 3 - Mr. Jeffrey R. Mannion

The requested information, or a request for an extension of time, should reference your above 510(k) number and should be submitted in duplicate to:

Food and Drug Administration Center for Devices and Radiological Health Document Mail Center (HFZ-401) 9200 Corporate Boulevard Rockville, Maryland 20850

If you have questions concerning the contents of this letter, please contact Donna Buckley at (301) 443-8243. If you need information or assistance concerning the IDE regulations, please contact the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its internet address "dsmo@fdadr.cdrh.fda.gov".

Sincerely yours,

Thomas J. Callahan, Ph.D.

Tou A kya for

Director

Division of Cardiovascular, Respiratory, and Neurological Devices Office of Device Evaluation Center for Devices and

Radiological Health



Mr. Jeffrey R. Mannion Regulatory Affairs Associate Angiodynamics, Incorporated 603 Queensbury Avenue Queensbury, New York 12804

Re: K964034

Transjugular Access Set Dated: September 30, 1996 Received: October 8, 1996

Dear Mr. Mannion:

We have reviewed your Section 510(k) notification of intent to market the device referenced above. We cannot determine if the device is substantially equivalent to a legally marketed predicate device based solely on the information you provided. To complete the review of your submission, we require the following additional information:

- 1. Please address the following issues regarding the indications for use:
  - a. Each indication for use needs to be specified and supported. You indicate in Section 5 that your device is "intended to be used for transjugular access in diagnostic and interventional procedures, such as biliary drainage, liver access and biopsy, placement of vena cava filters, and cholangiography." However, the device to which you claim substantial equivalence is indicated only for liver access. Therefore, please address the following: (1) revise the indication for your device to specify the areas where it is to be used; (2) identify a predicate device for each indication; and (3) demonstrate that your device is substantially equivalent to the predicate by providing comparative data and performance data.
  - b. Please list the vena cava filters that you intend to place with your device and provide either an explanation or testing to demonstrate the adequacy of your device to place these filters without incident.



## Page 2 - Mr. Jeffrey R. Mannion

- 2. The predicate device includes a pencil point needle, whereas your device includes a Ross needle and a beveled tip needle. Please provide either the trade names and manufacturers of devices that use the 14 gauge curved Ross needle and the 21 gauge beveled tip needle for the same indications or provide an article from a peer reviewed journal that indicates that these needle tips and gauges are adequate for the stated indications.
- 3. The package label indicates that the kit includes a 22 gauge beveled point needle, whereas the rest of the file references a 21 gauge beveled point needle. Please clarify this discrepancy.

We believe that this information is necessary for us to determine whether or not this device is substantially equivalent to a legally marketed predicate device with regard to its safety and effectiveness.

You may not market this device until you have provided adequate information described above and required by 21 CFR 807.87(k), and you have received a letter from FDA allowing you to do so. If you market the device without conforming to these requirements, you will be in violation of the Federal Food, Drug, and Cosmetic Act (Act). You may, however, distribute this device for investigational purposes to obtain clinical data if needed to establish substantial equivalence. Clinical investigations of this device must be conducted in accordance with the investigational device exemption (IDE) regulations.

If the information, or a request for an extension of time, is not received within 30 days, we will consider your premarket notification to be withdrawn and your submission will be deleted from our system. If you submit the requested information after 30 days it will be considered and processed as a new 510(k); therefore, all information previously submitted must be resubmitted so that your new 510(k) is complete.

Page 3 - Mr. Jeffrey R. Mannion

The requested information, or a request for an extension of time, should reference your above  $510\,(k)$  number and should be submitted in duplicate to:

Food and Drug Administration
Center for Devices and
Radiological Health
Document Mail Center (HFZ-401)
9200 Corporate Boulevard
Rockville, Maryland 20850

If you have questions concerning the contents of this letter, please contact Donna Buckley at (301) 443-8243. If you need information or assistance concerning the IDE regulations, please contact the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its internet address "dsmo@fdadr.cdrh.fda.gov".

Sincerely yours,

Thomas J. Callahan, Ph.D.
Director
Division of Cardiovascular, Respiratory,
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health



## Records processed under FOIA Request #2018-3544; Released by CDRH on 09-13-2018 \*\*DEPARTMENT OF HEALTH & HUMAN SERVICES\*\*

Page 4 - Mr. Jeffrey R. Mannion

Prepared by:DBUCKLEY:swf:1/2/97

final:swf:1/6/97 cc: HFZ-401 DMC

HFZ-404 510(k) Staff

HFZ-450 Division

D.O.

## FILE COPY

OFFICE	SURNAME	DATE	OFFICE	SURNAME	DATE	OFFICE	SURNAME	DATE
450	Buckley	1/6/97					1	
	Hamelon	16/97						
Ro	Kep	1/4/47						

U.S. GPO 1986-169-089




Revised:7-1-96

DEPARTMENT OF HEALTH & HUMAN SERVICES Public Health Service
Records processed under FOIA Request #2018-3544; Released by CDRH on 09-13-2018 Food And Drug Administration

J. C.		The second secon		Memorandum
From:	Reviewer(s) - Name(s) Donna B	uckley		
Subject:	510(k) Number	K96-	1034	
To:	The Record - It is my recommendation th	at the subject :	510(k) Notif	ication:
	☐ Is substantially equivalent to marketed	devices.		
	☐ Requires premarket approval. NOT sul	stantially equ	ivalent to m	arketed devices.
	☐ Requires more data.			
	☐ Accepted for review	·		
	(date)			
	☐ Other (e.g., exempt by regulation, not a	device, duplic	ate, etc.)	
Is this de	vice subject to Postmarket Surveillance?		□YES	⊠NO
Is this de	vice subject to the Tracking Regulation?		□YES	⊠NO
Was clini	cal data necessary to support the review of	this 510(k)?	□YES	<b>™</b> NO
Is this a p	rescription device?		<b>⊠</b> YES	□NO
Was this	510(k) reviewed by a Third Party?		□YES	⊠NO
Truthful (require	k) contains:  and Accurate Statement  Requested  Ered for originals received 3-14-95 and after)  c) summary OR A 510(k) statement quired certification and summary for class		Ā	:
☑ The inc	lication for use form (required for originals	received 1-1-	96 and after	)
The subm	itter requests under 21 CFR 807.95 (doesn	t apply for SE	is):	•
	nfidentiality   □ Confidentiality for 90 day		•	ntiality exceeding 90 days
Predicate	Product Code with panel and class: Add	itional Produc	t Code(s) wi	th panel (optional):
Review:_				
<b>(</b> B	ranch Chief) (B	ranch Code)		(Date)
Final Revi	ew:			
	(Division Director)			(Date)

## 510(k) "SUBSTANTIAL EQUIVALENCE" DECISION-MAKING PROCESS (DETAILED)



- SIO(k) submissions compare new devices to marketed devices. FDA requests additional information if the relationship between marketed and "predicate" (pre-Amendments or reclassified post-Amendments) devices is undear.
- \*\* This decision is normally based on descriptive information alone, but limited testing information is sometimes required.
- \*\*\* Data may be in the \$10(k), other \$10(k)s, the Center's classification files, or the literature.



## DEPARTMENT OF HEALTH & HUMAN SERVICES

**Public Health Service** Food and Drug Administration

## Memorandum

DATE:

January 2, 1997

TO:

File K964034

FROM:

Donna C. Buckley

Mechanical Engineer, FDA/CDRH/ODE/DCRND/ICDG

**SUBJECT:** 

510(k) Review of Transjugular Access Kit (catheter introducer) -- Angiodynamics

CONTACT: Brian Kunst -- (518) 798-1215 ext.190; FAX (518) 798-3625

DECISION: ADDITIONAL INFORMATION REQUIRED

INTENDED USE: The Transjugular Access Set is intended to be used for percutaneous transjugular access during diagnostic and/or therapeutic interventional procedures. The examples of procedures provided were biliary drainage, liver access and biopsy, placement of vena cava filters, and cholangiography. I spoke with Brenda Bolden (DDIG) and Richard Felton (DGRD) with regard to the device design and the indications for use. Both claimed that the indications for use need clarification and justification prior to a detailed review of the device design and testing. For example, there are no known predicate intravascular liver biopsy devices.

**DESCRIPTION:** The access kit is comprised of the following components and corresponding 510(k) numbers are noted where applicable:

- 1) 8F introducer sheath
- 2) hemostasis valve (K894446)
- 3) 8F PTFE catheter
- 4) 5F AngioOptic catheter (similar to K954423)
- 5) 14 guage curved Ross needle
- 6) 21 guage beveled tip needle
- 7) 1 cc syringe (K905447)
- 8) 10 cc syringe (K905447)
- 9) y-adapter (K905447)
- 10) dual check valve (K905447)

PREDICATE INFORMATION: The predicate device is Cook's Transjugular Liver Access Set. A letter was provided from Cook that stated that the device is a currently marketed modified preamendment device.

BENCH TESTING: Pull testing was performed on the 8F catheter, the 8F catheter introducer and the 5F AngioOptic angiographic catheter. Test results appear adequate.



**BIOCOMPATIBILITY:** Biocompatibility test information was provided for the following: (1) hemolysis; (2) cytotoxicity; (3) pyrogenicity; (4) irritation and (5) implantation. A systemic toxicity study was also performed with mice. Test results appear adequate.

## STERILIZATION/PACKAGING:

EtO sterilization will be performed in accordance with AAMI Guidelines for Industrial Ethylene Oxide Sterilization of Medical Devices (SAL > 10<sup>-6</sup>, ethylene oxide and ethylene chlorohydrin residuals < 25 ppm, ethylene glycol residuals < 250 ppm). Pyrogen testing will also be performed in accordance with the United States Pharmacopia Bacterial Endotoxins Test. The access set is packaged in Tyvek to polyethylene/polyester Mylar pouches. All packaging materials are identical to those used to package the AngioOptic angiographic catheters (K954423).

## **DOCUMENTATION:**

Indications for Use \_\_X\_

Truthful and Accurate Statement \_\_X\_

Safety and Effectiveness \_\_X\_ Statement \_\_\_\_ Summary

## **CORRESPONDENCE**:

Cook Regulatory Affairs -- April Lavender Ms. Lavander reported that their transjugular access kit is a modified preamendment device. Registration number: A176504.

## ADDITIONAL INFORMATION REQUESTED:

- 1. Please address the following issues regarding the indications for use:
  - a. Each indication for use needs to be specified and supported. You indicate in Section 5 that your device is "intended to be used for transjugular access in diagnostic and interventional procedures, such as biliary drainage. liver access and biopsy, placement of yena cava filters, and cholangiography."
    b. (b) (4)
    c. (b) (4)



	(b) (4)
2.	(b) (4)
	75) <i>(</i> 4)
3.	(b) (4)

## RECOMMENDATION:

I recommend that the above additional information is required. Several issues need to be addressed prior to a substantially equivalent determination.

- 1. The preamendment status of the predicate should be established.
- 2. A 510(k) kit certification should be obtained.
- 3. The indications for use should be clearly stated.
- 4. Depending on the indications, either consult requests or separate 510(k)s should be submitted to DRAERD and/or DDIG.
- 5. Testing should be evaluated for each indication.

Donva Butz

Donna Buckley

Consultation 16197

## DCRND Screening Checklist for Premarket Notification 510(k)

Device: Jansmonlas acress Cot	K969	4030
Submitter: Monodynamics	1 / 10	1001
Items which should be Included (circle missing & needed information)	Yes No	if Item Needed & MISSING
1. General information a) trade name, b) common name, c) establishment registration #, d) address of manufacturer, e) device class, f) new of modification, g) predicate device identified, h) 513/51/4 compliance (none yet available), i) Truth and Accuracy Statement, j) Indications for Use Statement (spparate page)	V	
SMDA requirements: 510(k) summary or statement (any Class device)	V	
Class III Certification & Summary (if Class III)	ALA	
3. Proposed Labeling: a) package labels, b) statement of intended use, c) advertisements or promotional materials d) MRI compatibility (if claimed)		
4. Description of device (or modification) including diagrams, engineering drawings, photographs, service manuals	V	
5. Comparison Information (similarities and differences) to named legally marketed equivalent device (table preferred) should include: a) labeling, b) intended use, c) physical characteristics, d) anatomical sites, f) performance (bench, animal, clinical) testing, g) safety characteristics		
<ol> <li>Biocompatibility data for all patient-contacting materials, OR, certification of identical material/formulation: a) component &amp; material, b) identify patient-contacting materials,</li> <li>c) biocompatibility of final sterilized product</li> </ol>		
<ol> <li>Sterilization and expiration dating information: a) sterilization method, b) SkL, c) packaging, d) specify pyrogen free,</li> <li>e) ETO residues, f) radiation dose</li> </ol>	1	· · · · · · · · · · · · · · · · · · ·
8. Software validation & verification: a) hazard analysis, b) level of concern, c) development documentation, d) certification	NA	
Meets current DCRND guidelines and applicable standards for this device: a) specify guidance, b) comply with content.		
Items shaded under "No" are necessary for all submissions. Circle checks in the "Needed & MISSING" column must be submitted be document.	cled items and it	tems with
ng: VesNo Reviewer: Show mallword	Date:/	0/11/96

## For DCRND Use Only

## DCRND Classification Checklist for Premarket Notification 510(k)

Device: Janspinglas acres Set Submitter: Angrodynamics	K964	034
Date received: Original 510(k): 8-47-96 This submission: 8-007-96	Revi cyc 1	le
Review Tier (circle one): I, II, III  (for Tier I, complete items 1-5 on the Screening Checklist)		
Question	Yes	No
A. Is the product a device?	X	
B. Is the device exempt from 510(k) by regulation or policy?		Х
C. Expedited Review Status: Requested by sponsor		X
Identified by DCRND		X
Granted by DCRND		
D. Has this device has been the subject of a previous NSE decision?		
If yes, does this new 510(k) address the NSE Issues(s), e.g., performance data?		
E. Has the sponsor been the subject of an integrity investigation?		
If yes, has the ODE Integrity Officer given permission to proceed with the review?		

Administrative Reviewer Signature

Date: 10/14/96



Food and Drug Administration Center for Devices and Radiological Health Office of Device Evaluation Document Mail Center (HFZ-401) 9200 Corporate Blvd. Rockville, Maryland 20850

October 08, 1996

ANGIODYNAMICS 603 QUEENSBURY AVE. QUEENSBURY, NY 12804 ATTN: BRIAN KUNST 510(k) Number: K964034 Received: 08-OCT-96

Product: TRANSJUGULAR ACCESS

SET

The Center for Devices and Radiological Health (CDRH), Office of Device Evaluation (ODE), has received the Premarket Notification you submitted in accordance with Section 510(k) of the Federal Food, Drug, and Cosmetic Act (Act) for the above referenced product. We have assigned your submission a unique 510(k) number that is cited above. Please refer prominently to this 510(k) number in any future correspondence that relates to this submission. We will notify you when the processing of your premarket notification has been completed or if any additional information is required. YOU MAY NOT PLACE THIS DEVICE INTO COMMERCIAL DISTRIBUTION UNTIL YOU RECEIVE A LETTER FROM FDA ALLOWING YOU TO DO SO.

On January 1, 1996, FDA began requiring that all 510(k) submitters provide on a separate page and clearly marked "Indication For Use" the indication for use of their device. If you have not included this information on a separate page in your submission, please complete the attached and amend your 510(k) as soon as possible. Also if you have not included your 510(k) Summary or 510(k) Statement, or your Truthful and Accurate Statement, please do so as soon as possible. There may be other regulations or requirements affecting your device such as Postmarket Surveillance (Section 522 (a) (1) of the Act) and the Device Tracking regulation (21 CFR Part 821). Please contact the Division of Small Manufacturers Assistance at the number below for more information.

Please remember that all correspondence concerning your submission MUST be sent to the Document Mail Center (HFZ-401) at the above letterhead address. Correspondence sent to any address other than the Document Mail Center will not be considered as part of your official premarket notification submission. Because of equipment and personnel limitations, we cannot accept telefaxed material as part of your official premarket notification submission, unless specifically requested of you by an FDA official. Any telefaxed material must be followed by a hard copy to the Document Mail Center (HFZ-401).

If you have procedural or policy questions, or want information on how to check on the status of your submission (after 90 days from the receipt date), please contact the Division of Small Manufacturers Assistance at (301) 443-6597 or their toll-free number (800) 638-2041, or call me at (301) 594-1190.

Sincerely yours,

Marjorie Shulman Consumer Safety Officer Premarket Notification Staff Office of Device Evaluation Center for Devices and Radiological Health



K964034

INCORPORATED

September 30, 1996

U.S. Food and Drug Administration Center for Devices and Radiological Health Document Mail Center (HFZ-401) 9200 Corporate Drive Rockville, Maryland 20850

Re:

510(k) Notification for

Transjugular Access Set

Dear Document Control Staff:

Pursuant to the requirements of Section 510(k) of the Food, Drug and Cosmetic Act, notification is made of the intention of AngioDynamics<sup>®</sup>, Inc. to manufacture and market the following:

## Transjugular Access Set

We consider our intent to market this device as confidential commercial information and request that it be considered as such by FDA. We have not disclosed our intent to market this device to anyone except employees of our establishment and have taken precautions to protect this confidentiality.

Thank you in advance for your consideration of this application. If there are any questions, please feel free to contact Brian Kunst, Director of RA/QA, at 518-798-1215, extension 123.

Sincerely,

ANGIODYNAMICS, Inc.

Jeffrey R. Mannion

Regulatory Affairs Associate

•		. PagcOI
510(k) Number (if known):		and complementations
Device Name: Transjugular Acce	ss Set	
·		· · · · · · · · · · · · · · · · · · ·
Indications For Use: The AngioDyna intended to b access during interventiona	e used 101   diagnostic	and/or therapeutic
		<b>,</b>
•		
(PLEASE DO NOT WRITE BELOW THI	S LINE - CONTINU	E ON ANOTHER PAGE IF NEEDED)
Concurrence of CDRH	, Office of Device	e Evaluation (ODE)
	Di ar	Division Sign-Off) ivision of Cardiovascular, Respiratory, and Neurological Devices K964634 IO(k) Number
Prescription Use	OR	Over-The-Counter Use
Per 21 CFR 801.109)		(Optional Format 1-2-96)
,		. (Optional connact -z->0)

H

## **TABLE OF CONTENTS**

<u>Contents</u>	<u>Section</u>	
Truthful And Accurate Statement	1	
General Information	2	
Summary and Certification	3	
Proposed labeling	4	
Device Description	5	
Comparative Information and Substantial Equivalence	6	
Biocompatability Assessment	7	
Sterilization Information	8	
ANGIODYNAMICS® Transjugular Access Set Drawings	9	
Predicate Device Labeling	10	



## **Truthful And Accurate Statement:**

[As required by 21 CFR 807.87(j)]

I certify that in my capacity as Regulatory Affairs Associate at AngioDynamics®, I believe to the best of my knowledge, that all data and information submitted in the premarket notification are truthful and accurate and that no material fact has been omitted.

Jeffrey R. Mannion

(Typed Name)

9/30/96

(Date)

(Premarket Notification (510(k) Number)



## General Information:

AngioDynamics, Inc. intends to introduce the following device into commercial distribution:

Transjugular Access Set Trade Name: a) Catheter Introducer **Device Name:** b) Established Registration Number: 319211 c) ANGIODYNAMICS® **Manufacturing Site Address:** d) 603 Queensbury Avenue Queensbury, New York 12804 Cosmed II, Inc. **Sterilization Site Address:** e) 8 Industrial Drive Coventry, Rhode Island 02816 1222273 Sterilizer Establishment f) Registration Number Classification: e)

870.1340 Catheter Introducer

Product Code: 74DYB

Class II 🔍 🐣

Cardiovascular Devices

f) **Device Equivalence:** 

This product is substantially equivalent to the following devices:

Cook® Rosch-Uchida Transjugular Liver Access Set

This device does not present additional risks to patients or different considerations regarding safety and effectiveness than those presented by the pre-Amendment devices.



g) Performance Standards: None Established

## h) ANGIODYNAMICS® Contact Information:

Name Brian Kunst, Director of RA/QA

Address: 603 Queensbury Avenue

Queensbury, New York 12804

Phone (518) 798-1215 extension 123

Fax: (518) 798-3625

## **Summary And Certification:**

I certify that in my capacity as Regulatory Affairs Associate at AngioDynamics<sup>®</sup>, I will make available all information included in this premarket notification on safety and effectiveness within 30 days of request by any person if the device described in the premarket submission is determined to be substantially equivalent. The information I agree to make available will be a duplicate of the premarket notification submission, including any adverse safety and effectiveness information, but excluding all patient identifiers and trade secrets and confidential commercial information as defined in 21 CFR 20.61.

Jeffrey R. Mannion

Regulatory Affairs Associate

Sol

**Proposed Labeling** 



## Instructions for Use



# PROPOSED ACCESS TECHNIQUE

8. Advance the 21 gauge beveled tip needle into the suspected location of the target site. Again, inject contrast to determine the exact locations of the 21

French catheter and the 21 gauge beveled needle, and advance the

5 French catheter through the 14 gauge Ross needle into the target

needle can be removed, being sure not to remove the Amplatz wire hen the proper device, depending on the procedure, can be placed nto the target site (i.e. biliary drainage catheter, balloon catheter, 9. Again, inject contrast to determine the exact locations of the 21 gauge needle and the target site.

30. When the correct location is verified, unlock the hubs of the 5

French catheter and the 21 gauge beveled needle, and advance the 5

French catheter through the 14 gauge Ross needle into the target site.

310.

311.

312.

314.

315.

316.

317.

318.

318.

319.

319.

319.

310.

310.

310.

310.

311.

311.

311.

311.

312.

313.

313.

314.

315.

315.

316.

317.

318.

318.

318.

319.

319.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310.

310. Once in the target site, the 5 French catheter and 14 gauge Ross

ANGIODYNAMICS\*

NCORPORATED

603 Queensbury Ave., Queensbury, N.Y. 12804 U.S.A. .el. 518-798-1215 FAX 518-798-3625

IC114 Rev. A 9/96

## **TRANSJUGULAR ACCESS SET**

## INSTRUCTIONS


FOR

USE

**ANGIODYNAMICS®** 

## SET ASSEMBLY

Figure 1 illustrates the proper set assembly.






Figure 2

# INDICATION FOR USE

- interventional diagnostic and/or therapeutic procedure. Percutaneous transjugular access to perform an
- The Access Set is supplied sterile and is intended for one time use only. Do not re-sterilize.

## CAUTIONS


- Federal law (U.S.A.) restricts this device to sale by or on he order of a physician.
- Store in a cool, dry place. Protect from UV light.
- aseptic techniques. Tighten all Luer lock fittings by hand Prior to use, assemble and disassemble the set using on the furnished catheters prior to using.

## WARNINGS

Transjugular access should be performed by a physician familiar with the technique, complications, and contraindications.



## COMPONENTS




14 Gauge Ross Needle



# 

		£ \$514,414,414,414,414,414,414,414,414,414,			######################################
indication for Use	Cautions		Components	Set Assembly	Proposed Access Technique
Ü	ĕ	W 33	Ö	Š	e E




## 

- sidearm of the hemostasis valve on the 8 French Introducer to Remove the 8 French catheter and inject contrast through the determine the exact location of the target site. 4
- Assemble the puncture complex (21 gauge beveled tip needle,



Ś



retracting the plunger of the 1cc syringe will draw physician to inject contrast in very small amounts to store the contrast, and use the Icc syringe as the injection syringe. By using the dual check valve, contrast from the 20cc syringe, thus enabling the

Package Label

## Transjugular Access Set

## Contents:

- 1-8F PTFE Introducer with Hemostasis Valve and Sidearm
- 1 8F PTFE Catheter
- 1 5F Nylon AngiOptic™ Catheter
- 1 14g Curved Ross Needle
- 1 22g Beveled Point Needle
- 1 1cc Syringe
- 1 10cc Syringe
- 1 Dual Check Valve
- 1 Hemostasis Y-Adapter

CATALOG NO. /Nº de catalogue /Número de catálogo /Kara:cg-Nr. /Numero di catalogo /Número do catálogo /Katalognummer /カタログ書号

QUANTITY /Quantité /Cantidad /Menge /Quantità /Quantidade :Antal /数量

LOT NO. /Nº de lot /Nº de lote /Posten-Nr. /Numero di lotto Lote nº /Lot nummer /ロット書号

STERILIZATION METHOD /Méthode de stérilisation .Metcoo de esterilización /Sterilisationsmethode /Metodo di sterilizzazione Metcoo de esterlização /Steriliseringsmetod /運置方法

STERILIZATION DATE /Date de stérilisation /Fecha de esterilización /Sterilisationsdatum /Data di sterilizzazione /Data de esterilização /Steriliseringsdatum /滅菌年月日

EXPIRATION DATE /Date de péremption Fecha de cacusidad Werfallsdatum /Data di scadenza /Data de expiração /Utgânçscatum 有効年月日

STERILE: Sterile unless package is opened or damaged. Single tatient use: Do not resterilize. CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

STERILE . SINGLE USE . STORE IN A COCL DRY PLACE . PROTECT FROM UV LIGHT Stérile · Destiné à un usage unique · Conserver dans un encroit sec et frais · Protéger de la lumière ultraviolette Esteril - Descrie a un usage unique - conserver dans che document rais - Proteger de la lura untravolette Esteril - Desse una sola vez - Guardar en un lugar sedo y fresco - Proteger de la luz ultravioleta Steril - Nur zum Einmaigeorauch - Kühl und trocken bewahren - Vor UV-Licht schützen Sterile - Uso unico solamente - Ricorre in un luggo fresco e secco - Protegere da la luce ultravioletto Estéril - Uso apenas uma vez - Deve ser conservado em che secco - Deve ser protegido contra radiação ultra violeta Steril - Endast for englaspociula - Forwas será con tort - Skyddas frân ultravolett fjus

滅菌済 再使用しないで下さい。 取扱説明書参照



## **Device Description:**

The AngioDynamics® Transjugular Access Set is intended to be used for transjugular access in diagnostic and interventional procedures, such as biliary drainage, liver access and biopsy, placement of vena cava filters, and cholangiography.

The AngioDynamics® Transjugular Access Set consists of the following components:

- 8 French introducer sheath with hemostasis valve
- 8 French PTFE catheter
- 5 French AngiOptic<sup>TM</sup> catheter
- 14 gauge curved Ross needle
- 21 gauge beveled tip needle




Fig. 5-1. Assembled Transjugular Access Set.

To allow for infusion of contrast agents to properly locate the device as it is advanced, the following components are also included:

• 1cc syringe, 10cc syringe, y-adapter, dual check valve

The introducer sheath consists of a 8 French I.D. x 12 French O.D. PTFE tube, 47cm in length with a nylon hub and a hemostasis valve. The intended use for the introducer sheath is to facilitate the placing of catheters into the vasculature, the same intended use as for all introducer sheaths. The introducer sheath is pictured in Fig. 5-2.




Fig. 5-2. Introducer Sheath with Hemostasis Valve.

63

The needle set consists of a 14 gauge blunt tipped 304 stainless steel curved Ross needle, and a 21 gauge beveled tip stainless steel needle with a nylon hub. The 14 gauge Ross needle has a torque handle / directional indicator on the hub. This handle is very useful for rotating the needle and to allow the user to know the direction in which the curve is pointing at all times. The 21 gauge needle is hollow to allow for contrast injections at any time to visualize location within the vasculature at any time throughout the procedure.

The intended use for the needle set is to act as a catheter stiffener, and to allow access to given areas within the vasculature. The needle set is pictured below in Fig. 5-3.




Fig. 5-3. 14 gauge curved Ross needle and 21 gauge needle with hub

The 8F catheter is comprised of a piece of PTFE tubing and a nylon catheter hub, and is intended to be used as a dilator with the Introducer Sheath (12F O.D., 8F I.D.) The 5F catheter is very similar to AngioDynamics® AngiOptic Catheter, K954423, and intended to be used to inject contrast into the patient during the procedure. The only difference between the AngiOptic Catheter used in the Transjugular Access Set and the straight

selective catheter which was found substantially equivalent in K954423 is the catheter length. The catheter in K954423 is 59.5cm in length while the 5 French straight selective catheter in K954423 is 90cm in length. The catheters are pictured in Fig. 5-4.





Fig. 5-4. 8F PTFE catheter and 5F AngiOptic catheter.







## Packaging:

The AngioDynamics® Transjugular Access Set is packaged in Tyvek to polyethylene\ polyester Mylar pouches. All outer packaging materials are identical to that used for AngioDynamics® AngiOptic Angiographic Catheter, 954423.



## Comparative Information and Substantial Equivalence:

The AngioDynamics® Transjugular Access Set is intended to be used for transjugular liver access in diagnostic and interventional procedures. This is the same intended use as the Cook Rösch-Uchida Transjugular Liver Access Set. Size and material comparison between the AngioDynamics® Transjugular Access Set and the Cook Rösch-Uchida Transjugular Liver Access Set is presented in Table 6-1.

	T	Cool Dre-1 Halid	
_	ANGIODYNAMICS	Cook Rösch-Uchida	
Features	Transjugular Access Set	Transjugular Liver	
		Access Set	
Components included	Introducer Sheath with	Introducer Sheath with	
	Hemostasis Valve	Hemostasis Valve	
	5F & 8F Catheter	5F & 10F Catheter	
	Curved stiffening needle	Curved stiffening needle	
	Access needle	Access needle	
	1cc Syringe		
	10cc Syringe		
	Dual Check Valve		
	Hemostasis Y-Adapter		
Introducer sheath inner diameter	8 French	10 French	
Introducer sheath material	PTFE	TFE	
Introducer sheath length	39.5cm	41.5cm	
R/O marker on sheath	No	Yes	
Catheter outer diameter	8 French, 5 French	10 French, 5 French	
Catheter lengths	8F - 51.5cm	10F - 51cm	
<u>-</u>	5F - 59.5cm	5F - 59cm	
Catheter material	PTFE	Radiopaque TFE	
Curved stiffening cannula size	14 gauge	14 gauge	
Curved stiffening cannula length	44.5cm	52cm	
Stiffening cannula material	Stainless steel	Stainless steel	
Access needle size/style	21 gage cannula	20 gage stylet	
-	beveled point	pencil point	
Access needle length	69cm	62cm	
Access needle material	Stainless steel	Stainless steel	

Table 6-1

The primary difference between the AngioDynamics® Transjugular Access System and the Cook Rösch-Uchida Transjugular Liver Access Set is the presence of the fluid management components (dual check valve, y-adapter, 1cc syringe, 10cc syringe), and 21 gauge hollow needle to permit contrast infusion for determination of needle location as it is being advanced. This will result in easier device location and will minimize the trauma created by inaccurate advancement.

The fluid management components have been found substantially equivalent in K905447, ANGIODYNAMICS® Pulse\*Spray Infusion Catheter System. The hemostasis valve was found substantially equivalent in K894446, Braun Hemostasis Valve. The 5 French AngiOptic™ catheter is very similar to the 5 French straight selective catheter, which was found substantially equivalent in ANGIODYNAMICS® AngiOptic Catheter, K954423. The only difference between the two catheters is the usable length. The 5 French straight selective catheter in K954423 is 90cm in length, while the 5 French AngiOptic™ catheter in the Transjugular Access Set is 59.5cm in length.

Pictorial comparison between the AngioDynamics<sup>®</sup> Transjugular Access Set and the Cook<sup>®</sup> Rosch-Uchida Transjugular Liver Access Set is given below in Fig. 6-1 and Fig. 6-2.




Fig. 6-1. AngioDynamics® Transjugular Access Set




Fig. 6-2. Cook® Rosch-Uchida Transjugular Liver Access Set

Enclosed as Attachment 6-1 is a letter from Cook®, Inc. stating that the Rosch-Uchida Transjugular Liver Access Set is a modified pre-Amendment device.



## Attachment 6-1



Records processed under FOIA Request #2018-3544; Released by CDRH on 09-13-2018

925 South Curry Pike P.O. Box 469

Bluenungton, IN 47402 U.S.A. Phone: 812 339 2236 Telex: 6711161 COOK UW Telefax: 812 339-8369





- - 82

June 21, 1995

Cook Incorporated

Ms. Maureeri Curtin 522 Woodland Trails East Stroudsburgh, PA 18301

**RE: REQUEST FOR 510(k) SUMMARIES** 

Dear Ms. Curtin:

This letter follows in response to your request for 510(k) Summaries of Safety and Effectiveness for the following items:

- Colapinto transjugular cholangiography and biopsy sets
- Ring transjugular intrahepatic access set
- Rösch-Uchida transjugular liver access set

These items are modified pre-Amendment devices, which are listed under the Listing Provisions of the Medical Device Amendments of 1976. This listing was required for those items that were already on the market prior to enactment of the Amendments in 1976, which called for premarket notifications. Consequently, no 510(k) summaries are available for these devices.

Thank you for this opportunity to be of service.

Sincerely,

COOK INCOMPORATED

April Lavender, RAC

Vice President


Regulatory Affairs

## Biocompatability assessment:

## Patient contact materials are:


- 304 stainless steel cannula tubing
- PTFE sheath/catheter tubing


PTFE is a very common material used in medical device applications and is well known for its biological inertness and lack of reactivity. The 304 stainless steel cannula tubing used in the 22 gauge beveled tip needle and the 14 gauge curved Ross needle is a very common material used in numerous medical devices throughout the industry. The 5 French AngiOptic Catheter biocompatibility test results from K954423 is included with all other biocompatibility test results in Attachment 7-1.





## Attachment 7-1


























































































Records processed under FOIA Request #2018-3544; Released by CDRH on 09-13-2018


	Pyrogen Test Specification
(b) (4)	



























































































PAGE \_\_\_\_1 OF \_\_10

DATE: 7/1/96

TO: Angio Dynamics

FAX: 518-798-3625


ATTN: Mark Steininger


FROM: Jeff Prioreschi


#### **COMMENTS:**


The information contained within this facaintile transmission is intended to be restricted to only the personal and confidential use of the addresses specified hereabove. This transmission is a confidential client communication and may also for other reasons be the addresses specified hereabove. This transmission is a confidential client communication and may also for other reasons be confidential or privileged. If you have received or otherwise obtained this transmitted information in error, then any review, reading, confidential or privileged. If you have received this transmission in error, please notify us immediately by telephone and return the complete original transmission to us by mail. Thank You.


ZEUS INDUSTRIAL PRODUCTS, INC.
501 BOULEVARD, SE, P.O. BOX 2167, ORANGEBURG, SOUTH CAROLINA 29116-2167
TOLL FREE: 800-526-3842 • PHONE: 803-531-2174 • FAX: 803-531-9244





























## Premarket Notification Transjugular Access Set

### **Sterilization Information:**



### Pyrogen Testing:


```
(b) (4)
```



## Premarket Notification Transjugular Access Set

**Product Drawings:** 






# Premarket Notification Transjugular Access Set

**Predicate Device Labeling:** 



#### RÖSCH-UCHIDA TRANSJUGULAR LIVER ACCESS SET

Used for transjugular liver access in diagnostic and interventional procedures. The radiopaque marker on the distal tip of the sheath enhances visualization during procedural use. Stiffening cannula provides added support during liver capsule puncture. Trocar stylet facilitates uncomplicated puncture. Supplied sterile in peel-open packages. Intended for one-time use.



SET ORDER Number	TROCAR STYLE		STIFFENING CANALLA			
		Diamete -	Length	Gage extra-thinwali	vergiti	
RUPS-100	, san aga	.038 inch (0.97 mni)	62.5 cm	14	51 to m	

REFERENCES: Page 283, R-BD32, R-BD3€

