

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS**REGULATORY AUTHORITY****JAN 0 8 2002**

Safe Medical Devices Act of 1990, 21 CFR 807.92

COMPANY NAME/CONTACT

Alan Curtis
 SURx
 2675 Collier Canyon Road
 Livermore, CA 94550
 (925) 398-4500 (phone)
 (925) 398-4509 (facsimile)
 acurtis@surx.com

NAME OF DEVICE

Trade Name:	<u>SURx LP System</u>
Common Name:	Electrosurgical System
Device Product Code:	GEI
Classification Name:	Electrosurgical Cutting and Coagulation Device and Accessories (21 CFR 878.4400)
Device Panel:	General Surgery/Restorative Devices
Device Classification:	Class II

PREDICATE DEVICES**Predicate Radiofrequency Electrosurgical Generator Systems**

- Vidamed TUNA[®] (Transurethral Needle Ablation) System (K960918, K965199, K002583)
- ArthroCare Urologic Multi-Electrode Electrosurgery System (K961069)
- Oratec Interventions Model ORA 40 Electrothermal Generator (K973159, K964071, K971532)

Predicate Indications for Use Devices

- Louisville Laboratories, Laparoscopic Bladder Neck Suspension Kit (K971797)
- Empi, Inc. Innova[®] ComfortPulse[®] Vaginal Electrode (K964577)

DEVICE DESCRIPTION

The SURx LP System consists of two components: the SURx LP Generator and the SURx LP Applicator. The SURx LP Applicator connects to the Generator. The Applicator also provides irrigation to the treatment site. The SURx LP Applicator is supplied sterile and intended for single use. The Applicator uses a bipolar design, which means that a return pad is not required for operation. The tip of the Applicator contains a thermistor to monitor temperature at the targeted area.

The SURx LP Generator is a radiofrequency (RF) electronic instrument.

Software is utilized in the operation of the SURx LP Generator.

INDICATION FOR USE STATEMENT

The SURx Laparoscopic Probe (LP) Radio Frequency (RF) System is indicated for shrinkage and stabilization of female pelvic tissue for treatment of Type II stress urinary incontinence due to hypermobility in women not eligible for major corrective surgery.

SUBSTANTIAL EQUIVALENCE COMPARISON

Technological Characteristics

The technological characteristics of the SURx LP System are identical to those of the cited predicate electrosurgical devices, as well as the devices distributed by other manufacturers. These devices are equivalent in terms of design, materials, principal of operation, and product specifications. Any differences between the SURx LP System and the predicate devices do not raise new issues regarding safety or effectiveness.

Indications for Use

Substantial equivalence is also supported for the SURx LP System by the predicate devices cleared for the treatment of female stress urinary incontinence.

Clinical Performance Data

Results of clinical evaluations were used to demonstrate that the SURx LP System functioned as clinically intended. Sufficient data have been gathered from clinical studies to determine that the SURx LP System performs as clinically intended and that no new issues of safety and effectiveness are introduced.

CONCLUSION

Based on the design, materials, function, intended use, and clinical evaluation, the SURx LP System is substantially equivalent to the devices currently marketed under the Federal Food, Drug and Cosmetic Act. In addition, the SURx LP System raises no new safety or effectiveness issues. Therefore, safety and effectiveness are reasonably assured, justifying 510(k) clearance for commercial sale.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. Alan Curtis
Vice President, Regulatory
and Clinical Affairs
SURx
2675 Collier Canyon Road
Livermore, California 94550

JAN 08 2002

Re: K011190

Trade/Device Name: SURx LP System

Regulation Number: 878.4400

Regulation Name: Electrosurgical cutting and coagulation device and accessories

Regulatory Class: II

Product Code: GEI

Dated: October 30, 2001

Received: October 31, 2001

Dear Mr. Curtis:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Mr. Alan Curtis

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read "Celia M. Witten". To the left of the signature is a small, stylized initial "CW".

Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

INDICATIONS FOR USE STATEMENT

510(k) Number: **K011190**

Device Name: SURx RF System

Indications for Use:

The SURx Laparoscopic Probe (LP) Radio Frequency (RF) System is indicated for shrinkage and stabilization of female pelvic tissue for treatment of Type II stress urinary incontinence due to hypermobility in women not eligible for major corrective surgery.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21 CFR 801.109)

or Over-The-Counter Use _____

Steph Rhoads
(Division Chief)
Division of General, Restorative
and Neurological Devices

510(k) N K011190