

JAN 16 2002

K013833

510(k) SUMMARY

The 510(k) Summary is submitted in accordance with 21 CFR §807.92 and the requirements of the SMDA of 1990.

1. Submitter's Name: Guidant Corporation
Vascular Intervention
2. Submitter's Address: 26531 Ynez Road
3. Telephone: (909) 914-4581
4. Fax: (909) 914-0339
5. Contact Person: Jennifer Pae Riggs
6. Date Prepared: November 16, 2001
7. Device Trade Name: HI-TORQUE BALANCE MIDDLEWEIGHT™
UNIVERSAL Guide Wires with Hydrocoat
Hydrophilic Coating
8. Device Common Name: Guide Wire
9. Device Classification Name: Catheter Guide Wire (74DQX)
10. Predicate Device: HI-TORQUE BALANCE MIDDLEWEIGHT™ Guide
Wire with Hydrocoat Hydrophilic Coating (K973494)
HI-TORQUE WHISPER™ MS Guide Wire with
Hydrocoat Hydrophilic Coating (K002206)
11. Device Description:

The HI-TORQUE BALANCE MIDDLEWEIGHT™ UNIVERSAL Guide Wire is a steerable guide wire available in a maximum diameter of 0.0137" and in lengths of 175 cm, 190 cm and 300 cm. The distal segment of the guide wire, up to the hypotube, is coated with hydrocoat to reduce friction for improved guide wire movement within the catheter. The distal tip is offered in a straight shapeable configuration and a pre-shaped "J" configuration. The proximal core has a maximum diameter of 0.0137". The proximal end of the guide wire is coated with PTFE, which reduces friction of the wire within a catheter. The BMW™ Universal Guide Wire is DOC® extendable in the 175 cm and 190 cm lengths.

12. Intended Use:

To facilitate the placement of balloon dilatation catheters during percutaneous transluminal coronary angioplasty (PTCA) and percutaneous transluminal angioplasty (PTA). The wire is also intended to facilitate the placement of compatible stent devices during therapeutic intravascular procedures.

13. Technological Characteristics:

Comparisons of the new and predicate devices show that the technological characteristics such as product performance, design and intended use are substantially equivalent to the current marketed predicate devices.

14. Performance Data:

In vitro bench testing and *in vivo* performance evaluations demonstrated that the HI-TORQUE BALANCE MIDDLEWEIGHT™ UNIVERSAL Guide Wire met the acceptance criteria and performed similarly to the predicate devices. No new safety or effectiveness issues were raised during the testing program and therefore, the HI-TORQUE BALANCE MIDDLEWEIGHT™ UNIVERSAL Guide Wire may be considered substantially equivalent to the predicate devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Ms. Jennifer Pae Riggs
Sr. Regulatory Affairs Coordinator
Guidant Corporation
26531 Ynez Road
Temecula, CA 92591-4628

JAN 16 2002

Re: K013833
HI-TORQUE BALANCE MIDDLEWEIGHT™ UNIVERSAL Guide Wires with
Hydrocoat Hydrophilic Coating
Regulation Number: 870.1330
Regulation Name: catheter guide wire
Regulatory Class: Class II
Product Code: DQX
Dated: November 16, 2001
Received: November 19, 2001

Dear Ms. Riggs:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality


Page 2 - Ms. Jennifer Pae Riggs

systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4586. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



Bram Zuckerman, M.D.
Acting Director
Division of Cardiovascular and
Respiratory Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

INDICATION FOR USE STATEMENT

510(k) Number (if known): K013833

Device Names: HI-TORQUE BALANCE MIDDLEWEIGHT™ UNIVERSAL Guide Wire with Hydrocoat Hydrophilic Coating

Indications for Use: To facilitate the placement of balloon dilatation catheters during percutaneous transluminal coronary angioplasty (PTCA), percutaneous transluminal angioplasty (PTA) and compatible stent devices during therapeutic intravascular procedures.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use OR Over-The-Counter
(Per 21 CFR 801.109) (Optional Format 1-1-96)

[Signature]
Division of Cardiovascular & Respiratory Devices
510(k) K013833