

JUN 25 1996

K 961390

Appendices

510(k) Summary of Safety and Effectiveness (App. A)

Statement Information supporting claims of substantial equivalence, as defined under the Federal Food, Drug and Cosmetic Act, respecting safety and effectiveness is summarized below. For the convenience of the Reviewer, this summary is formatted in accordance with the Agency's final rule "...510(k) Summaries and 510(k) Statements..." (21 CFR §807) and can be used to provide a substantial equivalence summary to anyone requesting it from the Agency.

Device description The ENDOPATH® ETS Endoscopic Linear Cutter and ETS-FLEX Endoscopic Articulating Linear Cutter devices deliver two double-staggered or two triple-staggered rows of staples while simultaneously dividing the tissue between the rows. A "no-knife" version may also be provided. The instrument's safety lock-out feature is designed to prevent a spent reload from being refired. Both the standard and the vascular/thin instruments have a staple line that is 20-60mm long and a cut line of 18-58mm long. A staple retaining cap on the reload protects the staple leg points during shipping and transportation. An articulation lever on the ETS-FLEX Endoscopic Articulating Linear Cutter enables bilateral movement of the instrument jaws.

The instrument is reloadable with either a standard, blue reload for tissue that is compressible to 1.5mm in thickness, a thick, green reload for tissue that is compressible to 2.0mm in thickness, or a vascular/thin, white reload for tissue that is compressible to 1.0mm in thickness. Do not reload the instrument more than seven times for a maximum of eight firings per instrument.

Intended use For transection, resection, and/or creation of anastomoses.

Indications statement The ENDOPATH® ETS Endoscopic Linear Cutter and ETS-FLEX Endoscopic Articulating Linear Cutter have application in general, urologic, gynecologic, and thoracic surgery for transection, resection, and/or creation of anastomoses.

Continued on next page

510(k) Summary of Safety and Effectiveness (App. A), Continued

Technological characteristics The technological characteristics of the New Devices are the same as the Predicate Device.

Performance data Pre-clinical laboratory evaluations were performed to ensure that the device can be used as designed. The studies demonstrated acceptable performance in transecting, resecting, and/or creation of anastomoses.

Conclusion **Based on the 510(k) summaries and 510(k) statements (21 CFR §807) and the information provided herein, we conclude that the New Devices are substantially equivalent to the Predicate Device under the Federal Food, Drug and Cosmetic Act.**

Contact Lorri Chavez
Ethicon Endo-Surgery, Inc.
4545 Creek Road
Cincinnati, Ohio 45242

Date April 9, 1996
