

K961488



Diagnostics

JUN 12 1996

510(k) Summary

Introduction According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.

1. Submitter name, address, contact Boehringer Mannheim Corporation
2400 Bisso Lane
P.O. Box 4117
Concord, CA 94520-4117
(510) 674 - 0690, extension 8415

Contact Person: Mary Koning

Date Prepared: April 9, 1996

2. Device name Proprietary name: Elecsys® T-Uptake Assay
Common name: Electrochemiluminescence assay for the thyroxine-binding capacity (T-Uptake).

Classification name: Thyroxine-binding globulin test system

3. Predicate device We claim substantial equivalence to the Enzymun-Test® TBK (T4 Uptake)

4. Device Description The Elecsys® test principle is based on competition principle. Total duration of assay: 18 minutes (37° C).
•1st incubation (9 minutes): Sample (15 µL), exogenous T4, and biotinylated T4-polyhapten (75 µL).
•2nd incubation (9 minutes): After addition of a specific anti-T4 antibody labeled with a ruthenium complex (75 µL), the polyhapten and the antibody derivative react to form a complex, the concentration of which is inversely proportional to the concentration of the excess, exogenous T4. This immunological complex is bound to the added streptavidin-coated microparticles (35 µL) via the interaction of biotin and streptavidin.

Continued on next page



510(k) Summary, Continued

4.
Device
Description,
cont.

•The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier (0.4 second read frame).

•Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the reagent bar code.

5.
Intended use

Immunoassay for the in vitro quantitative determination of thyroxine-binding capacity in human serum and plasma.

Continued on next page

510(k) Summary, Continued

6. Comparison to predicate device

The Boehringer Mannheim Elecsys® T-Uptake Assay is substantially equivalent to other products in commercial distribution intended for similar use. Most notably it is substantially equivalent to the currently marketed Enzymun-Test® TBK (T4 Uptake).

The following table compares the Elecsys® T-Uptake with the predicate device, Enzymun-Test® TBK (T4 Uptake). Specific data on the performance of the test have been incorporated into the draft labeling in attachment 5. Labeling for the predicate device is provided in attachment 6.

Similarities:

- Detection of Thyroxine-Binding Capacity (T-Uptake)
- Sample type: Serum and plasma
- Reportable range of 0.2 - 1.9 TBI
- Standardized according to Enzymun-Test® TBK Assay
- Polyclonal Antibody: Sheep polyclonal anti-T-Uptake
- Solid phase binding principle: Streptavidin/Biotin

Differences:

Feature	Elecsys® T Uptake	Enzymun-Test® TBK
Reaction test principle	Electrochemiluminescence	ELISA/1-step sandwich assay using streptavidin technology
Instrument required	Elecsys® 2010	ES 300
Calibration Stability	A calibration is recommended every 7 days if kit is not consumed; 4 weeks with same reagent lot if reagent is consumed within 7 days.	Calibration required every run

Continued on next page

510(k) Summary, Continued

6.
Comparison
to predicate
device, (cont.)

Performance Characteristics:

Feature	Elecsys® T Uptake	Enzymun-Test® TBK
Precision	Modified NCCLS (TBI):	Modified NCCLS (TBI):
Level	Sample Control 1 Control 2	Low Mid High
N	60 60 60	118 120 119
Within-Run	0.95 0.94 1.09	0.61 0.93 1.28
%CV	2.15 3.29 2.39	5.7 2.6 1.4
Total	0.95 0.94 1.09	0.61 0.93 1.28
%CV	3.25 3.68 2.66	7.2 3.4 2.0
Lower Detection Limit	0.2 TBI	0.2 TBI
Linearity	0.2 TBI to the value of the highest standard	0.2 TBI to the value of the highest standard
Method Comparison	Vs Enzymun-Test® TBK <u>Least Squares</u> $y = 0.99x - 0.03$ $r = 0.908$ $SEE = 0.04$ $N = 319$ <u>Passing/Bablok</u> $y = 1.07x - 0.12$ $r = 0.908$ $SEE = 0.04$ $N = 319$	Vs Enzymun-Test® TBK <u>Least Squares</u> $y = 1.02x + 0.084$ $r = 0.986$ $SEE = 0.054$ $N = 52$

Continued on next page



510(k) Summary, Continued

6. Comparison to predicate device, (cont.)

Performance Characteristics, cont.:

Feature	Elecsys® T Uptake	Enzymun-Test® TBK
Interfering substances	No interference at:	No interference at:
Bilirubin	25 mg/dL	64.5 mg/dL
Hemoglobin	1 g/dL	1 g/dL
Lipemia	1500 mg/dL	1250 mg/dL
Biotin	30 ng/mL	200 ng/mL
Specificity	% Cross-reactivity	% Cross-reactivity
L-T4	100	100
D-T4	100	100
L-T3	1.5	3.5
D-T3	1.4	2.9
3-iodo-L-tyrosine	0.002	<0.1
3,5-diiodo-L-tyrosine	0.01	<0.1
Tetraiodo-thyroacetic acid	38.5	20