

K964068  
Dec. 13, 1996

### SUMMARY OF SAFETY AND EFFECTIVENESS

#### Carbamazepine Method for Bayer Technicon Immuno 1<sup>®</sup> System

Listed below is a comparison of the performance between the Immuno 1 Carbamazepine method and a similar device that was granted clearance of substantial equivalence (Syva EMIT<sup>®</sup> 2000 Carbamazepine Assay, Behring Diagnostics Inc.). The information used in the Summary of Safety and Effectiveness was extracted from the Immuno 1 Carbamazepine method sheet and the Syva EMIT<sup>®</sup> 2000 Carbamazepine Assay insert.

#### INTENDED USED

This *in vitro* method is intended to quantitatively measure carbamazepine, an anti-convulsant drug, in human serum or plasma (heparin) using Syva EMIT<sup>®</sup> 2000 Assay on a *Technicon Immuno-1<sup>®</sup>* system. Measurements of carbamazepine are used in the diagnosis and treatment of carbamazepine overdose and in monitoring serum levels of carbamazepine to ensure appropriate therapy.

METHOD	Immuno 1 Carbamazepine	Syva EMIT <sup>®</sup> 2000 (predicate Device)
Part No.	Reagents T01-3677-51 Calibrators T03-3679-01	Reagents 4F019UL Calibrators 4F109UL
Minimum Detectable Conc.	0.2 µg/mL	0.5 µg/mL
Precision (Between-Run)	2.8 µg/mL 3.9% 9.6 µg/mL 4.2% 16.1 µg/mL 8.7%	4.0 µg/mL 7.3% 10.0 µg/mL 5.9% 16.3 µg/mL 5.6%
Correlation	$y = 1.07x - 0.13$	

where

y = Immuno 1 Carbamazepine method

x = Syva EMIT<sup>®</sup> 2000 Carbamazepine Assay\*

n = 50

r = 0.978

S<sub>yx</sub> = 0.89 µg/mL

\*This assay was performed on Roche COBAS FARA II<sup>®</sup> Instrument.



Gabriel J. Murasa, Jr.  
Manager Regulatory Affairs  
Bayer Corporation  
511 Benedict Avenue  
Tarrytown, New York 10591-5097

  
Date