

**Boston
Scientific
Corporation**

One Boston Scientific Place
Natick, MA 01760-1537
508-650-8000

K964284
JAN 4, 1997

510(k) Summary

Pursuant to §513(l)(3)(A) of the Food, Drug, and Cosmetic Act, Boston Scientific Corporation is required to submit with this Premarket Notification "... adequate summary of any information respecting safety and effectiveness or state that such information will be made available upon request of any person." Boston Scientific Corporation chooses to submit a summary of information respecting safety and effectiveness.

Trade Name: Stainless Steel Greenfield®
Vena Cava Filter with
12F Introducer System

Manufacturer: Boston Scientific Corporation
480 Pleasant Street
Watertown, MA 02172

Device Generic Name: Guidewire for Vena Cava Filter System

Classification:

According to Section 513 of the Federal Food, Drug and Cosmetic Act, FDA has classified cardiovascular intravascular filters in Class III, 21 CFR § 870.3375, under guidance by the Circulatory Systems Devices Panel. Although 21 CFR § 870.3375 requires premarket approval for cardiovascular intravascular filters, no effective date has been established for the requirement of PMA's for these devices at this time.

Predicate Device: Stainless Steel Greenfield®
Vena Cava Filter with
12F Introducer System

Device Description:

The **Stainless Steel Greenfield® Vena Cava Filter with 12F Introducer System** consists of sequential vascular dilators for percutaneous access; a guidewire to guide the Introducer Sheath and Dilator into place; an Introducer Sheath and Dilator to create a passageway to guide the Introducer Catheter to the implant site; and an Introducer Catheter to transport the pre-loaded filter through the sheath and over the guidewire to the implant site for filter deployment .

Indications For Use:

The **Stainless Steel Greenfield® Vena Cava Filter with 12F Introducer System** is indicated for the prevention of pulmonary embolism via placement in the vena cava in the following situations:

68

- Venous thrombosis or pulmonary thromboembolism when anticoagulants are contraindicated or inadequate for management of venous thrombosis with significant risk of, or following, pulmonary thromboembolism.
- Failure of anticoagulant therapy in thromboembolic disease.
- Emergency treatment following massive pulmonary embolism where anticipated benefits of conventional therapy are reduced.
- Chronic, recurrent pulmonary embolism where anticoagulant therapy has failed or is contraindicated.

Contraindications:

Patients in whom the diameter of the inferior vena cava exceeds 28 mm are contraindicated for the placement of the **Stainless Steel Greenfield® Vena Cava Filter with 12F Introducer System**.

- Presence of thrombus at the femoral puncture site, in the iliac vein, or in the inferior vena cava risks the dislodgement of thrombus during catheter manipulation.
- Patients in whom pregnancy has been confirmed.

Safety and Performance:

Bench and Animal testing were performed and the data supported the substantial equivalence of the proposed guidewire for the **Jugular Stainless Steel Greenfield® Vena Cava Filter with 12F Introducer System** to the current guidewire of the **Jugular Stainless Steel Greenfield® Vena Cava Filter with 12F Introducer System**.

The biocompatibility of the **Stainless Steel Greenfield® Vena Cava Filter with 12F Introducer System (K912038)** was previously established through biocompatibility testing performed on the predicate in accordance with the Tripartite Biocompatibility Guidance (#G87-1). This testing is also in accordance with blue book memorandum #G95-1, "Use of International Standard ISO-10993, Biological Evaluation of Medical Devices Part 1: Evaluation and Testing with FDA-modified matrix."

Conclusion:

Based upon the indications for use, technological characteristics, and safety and performance testing, the guidewire for the **Jugular Stainless Steel Greenfield® Vena Cava Filter with 12F Introducer System** meets the requirements for its intended use.

169