

K964442

DEC 18 1996

SUMMARY OF SAFETY AND EFFECTIVENESS**Valproic Acid Method for Bayer Technicon Immuno 1<sup>®</sup> System**

Listed below is a comparison of the performance between the Immuno 1 Valproic Acid method and a similar device that was granted clearance of substantial equivalence (Syva EMIT<sup>®</sup> Valproic Acid Assay, Behring Diagnostics Inc.). The information used in the Summary of Safety and Effectiveness was extracted from the Immuno 1 Valproic Acid method sheet and the Syva EMIT<sup>®</sup> Valproic Acid Assay Insert Sheet.

**INTENDED USED**

This *in vitro* method is intended to quantitatively measure valproic acid, an anticonvulsant drug, in human serum or plasma (heparin) using Syva EMIT<sup>®</sup> Valproic Acid Assay on a *Technicon Immuno-1<sup>®</sup>* system. Measurements of valproic acid are used in the diagnosis and treatment of valproic acid overdose and in monitoring levels of valproic acid to ensure appropriate therapy.

METHOD	Immuno 1 Valproic Acid Method Set	Syva EMIT <sup>®</sup> Valproic Acid Assay (predicate Device)
Part No.	T01-3698-01	6G009UL
Minimum Detectable Conc.	0.9 µg/mL	Not listed
Precision (Between-Run)	3.4% @ 50.4 µg/mL 2.6% @ 115.4 µg/mL 2.9% @ 132.2 µg/mL	4.6% @ 76 µg/mL
Correlation	$y = 0.94x - 0.22$  where $y =$ Immuno 1 Valproic Acid method $x =$ Syva EMIT <sup>®</sup> Valproic Acid Assay* $n =$ 50 $r =$ 0.996 $S_{yx} =$ 3.16 µg/mL	

\*This assay was performed on COBAS FARA II<sup>®</sup> Instrument using parameters and protocol specified in Behring Application Sheet.

Gabriel J. Muraca, Jr.  
 Manager Regulatory Affairs  
 Bayer Corporation  
 511 Benedict Avenue  
 Tarrytown, New York 10591-5097

11/5/96  
 Date