

**510(k) SUMMARY**  
**Summary of Safety and Effectiveness**

**APPLICANTS NAME AND ADDRESS:**

Drager Inc.  
Critical Care Systems  
4101-100 Pleasant Valley Road  
Chantilly, VA 20151

AUG 26 1997

**APPLICANTS TELEPHONE NUMBER:**

(703)-817-0100

**APPLICANTS FACSIMILE NUMBER:**

(703)-817-0101

**APPLICANTS CONTACT PERSON:**

Harald Kneuer  
Regulatory Affairs Manager

**DATE THE SUMMARY WAS PREPARED:**

November 11, 1996

**DEVICE NAME:**

Trade Name:	Evita 2 dura
Common Name:	Evita 2 dura
Classification Name:	Ventilator, Continuous (per 21 CFR 868.5895)

**LEGALLY MARKETED DEVICE TO WHICH DRAGER INC. IS  
CLAIMING SUBSTANTIAL EQUIVALENCE:**

**Evita 4** - Manufactured by Dragerwerk AG, Lubeck, Germany and sold in the United States by Drager Inc.

**Evita** - Manufactured by Dragerwerk AG, Lubeck, Germany and sold in the United States by Drager Inc.

**Servo 300** - Sold in the United States by Siemens

## Intensive Care Ventilator Data Comparison Sheet

Device	Evita 4	Evita 2 dura	Evita	Servo 300
<b>510 (K) number</b>	K 961687		K 924930	K902859
<b>Manufacturer</b>	Dräger	Dräger	Dräger	Siemens
<b>Classification</b>	Intensive Care Ventilator	Intensive Care Ventilator	Intensive Care Ventilator	Intensive Care Ventilator
<b>Working Principle</b>	time cycled, volume constant, pressure controlled.	like Evita 4 like Evita 4	time cycled, volume constant, pressure controlled.	time cycled, volume constant, pressure controlled.
<b>Software controlled</b>	Yes	Yes	Yes	Yes
<b>Gas delivery system</b>	Fast "High Pressure Servo Valve", microprocessor controlled	like Evita 4	Fast "High Pressure Servo Valve", microprocessor controlled	Fast "High Pressure Servo Valve", microprocessor controlled
<b>Voltage for gas delivery system</b>	12 V DC	like Evita 4	12 V DC	-

Device	Evita 4	Evita 2 dura	Evita	Servo 300
Pneumatic Unit	Partnumber 84 12 360	like Evita 4		
Software for Pneumatic Unit	Partnumber 84 11 742	like Evita 4		
Electronic Unit	Partnumber 84 12 835	like Evita 4		
Software for Electronic Unit	Partnumber 84 11 741	like Evita 4		
<b>Modes of Ventilation:</b>				
CMV	Yes	Yes	Yes	Yes
SIMV	Yes	Yes	Yes	Yes
CPAP	Yes	Yes	Yes	Yes
ASB	Yes	Yes	Yes	Yes
Assist Control	Yes	Yes	Yes	Yes
MMV	Yes	Yes	Yes	No
PCV+ (BIPAP)	Yes*1	Yes*1	Yes	Yes
APRV	Yes	Yes	Yes	No
CMV, Assist Control, SIMV and MMV with AutoFlow	Yes*2	Yes*2	No	Yes
ILV	Yes	Yes	Yes	

\* see explanation of differences

Device	Evita 4	Evita 2 dura	Evita	Servo 300
<b>Ventilator Settings</b>				
Rate	0.5 to 100 bpm	like Evita 4	5 to 60 bpm for CMV 0.5 to 20 bpm for SIMV	5 to 150 bpm for CMV 0.5 to 40 bpm for SIMV
Vt	0.02 to 0.3 L (Pediatric Mode)* <sup>3</sup> 0.1 to 2.0 L (Adult Mode)* <sup>3</sup>	like Evita 4	0.1 to 2.0 L	0.002 to 4.0 L
T <sub>insp</sub>	0.1 ... 10 s	like Evita 4	not directly adjustable	10 to 80 % (stepless)
I:E-Ratio	not directly adjustable* <sup>4</sup>	like Evita 4	4:1 to 1:5	not directly adjustable
Flow	6 to 30 lpm (Paediatric Mode)* <sup>3</sup> 6 to 120 lpm (Adult Mode)* <sup>3</sup>	like Evita 4	6 to 120 lpm	not directly adjustable
P <sub>max</sub>	0 to 100 cmH <sub>2</sub> O	like Evita 4	0 to 100 cmH <sub>2</sub> O	10 to 120 cmH <sub>2</sub> O
PEEP	0 to 35 cmH <sub>2</sub> O	like Evita 4	0 to 35 cmH <sub>2</sub> O	0 to 50 cmH <sub>2</sub> O
O <sub>2</sub> -Conc.	21 to 100 Vol. %	like Evita 4	21 to 100 Vol. %	21 to 100 Vol. %
<b>Flow Trigger Sensivity</b>				
for CMV assist	0.3 to 15 lpm* <sup>5</sup>	like Evita 4	0.5 to 5.0 cmH <sub>2</sub> O	0 to 17 cmH <sub>2</sub> O P. trigger
for SIMV	0.3 to 15 lpm* <sup>5</sup>	like Evita 4	0.7 cmH <sub>2</sub> O fixed trigger	0 to 17 cmH <sub>2</sub> O P. trigger
for ASB (pressure support)	0.3 to 15 lpm	like Evita 4	1 to 15 lpm	0 to 17 cmH <sub>2</sub> O P. trigger
for CPAP	0.2 cmH <sub>2</sub> O fixed trigger	like Evita 4	0.2 cmH <sub>2</sub> O fixed trigger	
P <sub>insp</sub>	0 to 80 cmH <sub>2</sub> O	like Evita 4	10 to 100 cmH <sub>2</sub> O	0 to 100 cmH <sub>2</sub> O
Rise time	0 to 2 s	like Evita 4	0 to 2 s	0 to 10 %
Pressure Support	0 to 80 cmH <sub>2</sub> O	like Evita 4	0 to 80 cmH <sub>2</sub> O	0 to 100 cmH <sub>2</sub> O

\* see explanation of differences

Device	Evita 4	Evita 2 dura	Evita	Servo 300
<b>Parameter setting</b>				
Main ventilation settings	with 10.4" Color TFT-LC-display, Touch screen and central	With keys, 7-Seg. Display rotary knob and central rotary dial knob for selecting and confirming* <sup>5</sup>	rotary knob for every parameter	rotary knob for every parameter
further settings	With 10.4" Color TFT-LC-Display, Touch screen and central rotary knob	Via menu	menu via hardkeys	-
<b>Monitors/Indicators</b>				
Power ON	TFT-LC-Display (Color)	LC-Display (monochrom) or TFT (Color)	LC-Display	LED
Dimension of Display	21.1cm x 15.8cm 640 x 480 dots	12.5cm x 9.5cm 640 x 480 dots	9.8cm x 7.9cm 128 x 240 dots	
Alarm Condition	red LED and TFT- LC-Display	red LED and TFT or LC-Display	red LED and LC Display	red LED and LED Display
Warnings	yellow LED and TFT-LC-Display	yellow LED and LC or TFT	yellow LED and LC-Display	
Messages	TFT-LC-Display	LCD-or TFT-Display	LC-Display	LED-Display
Curves	TFT-LC-Display	LCD-or TFT-Display	LC-Display	no

\* see explanation of differences

Device	Evita 4	Evita 2 dura	Evita	Servo 300
Measured values	TFT-LC-Display	LC-or TFT Display	LC-Display and 7-Segment-Display	several LED-Displays
Adjusted values	TFT-LC-Display	7-Segment Display and LC- or TFT- Display	Scale of potentiometer and LC-Display	Scale of potentiometer and LED-Display
Audible Alarm	Yes	Yes	Yes	Yes
Alarm sound level	adjustable	adjustable	fixed	adjustable
<b>Pressure Limitation</b>				
adjustable Pmax	0 to 100 cmH <sub>2</sub> O	like Evita 4	0 to 100 cmH <sub>2</sub> O	10 to 100 cmH <sub>2</sub> O
Safety Valve	100 cmH <sub>2</sub> O	like Evita 4	100 cmH <sub>2</sub> O	120 cmH <sub>2</sub> O
Inspiratory relief valve	Yes	like Evita 4	Yes	Yes

Device	Evita 4	Evita 2 dura	Evita	Servo 300
<b>Alarms</b>				
Airway Pressure High /Low	Yes	Yes	Yes	Yes
Minute volume High / Low	Yes	Yes	Yes	Yes
Apnea	Yes	Yes	Yes	Yes
Gas Supply Failure	Yes	Yes	Yes	Yes
Power failure	Yes	Yes	Yes	Yes
Alarm Silence / Reset	Yes	Yes	Yes	Yes
Ventilator Inoperative	Yes	Yes	Yes	Yes
Breathing frequency High	Yes	Yes	Yes	No
Oxygen concentration High / Low	Yes	Yes	Yes	Yes
Breathing gas temperatur High	Yes	Yes	Yes	No
Endexpiratory CO <sub>2</sub> - concentration etCO <sub>2</sub> High / Low	Yes	Yes* <sup>8</sup> (option)	No	No
O <sub>2</sub> saturation SpO <sub>2</sub> High / Low	Yes (option)* <sup>9</sup>	Yes* <sup>9</sup>	No	No
Pulse High / Low	Yes (option) * <sup>9</sup>	Yes* <sup>9</sup>	No	No

\* see explanation of differences

Device	Evita	Evita 2 dura	Evita	Servo 300
<b>Displayed Values</b>				
max. airway pressure P <sub>peak</sub>	0 to 99 cmH <sub>2</sub> O	like Evita 4	0 to 99 cmH <sub>2</sub> O	Yes
Plateau pressure	0 to 99 cmH <sub>2</sub> O	like Evita 4	0 to 99 cmH <sub>2</sub> O	Yes
PEEP	0 to 99 cmH <sub>2</sub> O	like Evita 4	0 to 99 cmH <sub>2</sub> O	Yes
mean airway pressure	0 to 99 cmH <sub>2</sub> O	like Evita 4	0 to 99 cmH <sub>2</sub> O	Yes
min. airway pressure	0 to 99 cmH <sub>2</sub> O	like Evita 4	0 to 99 cmH <sub>2</sub> O	No
insp. oxygen concentration	15 to 100 Vol.-%	like Evita 4	5 to 99 Vol.-%	Yes
Minute volume MV	0 to 99 lpm	like Evita 4	0 to 99 lpm	Yes
spont. breath. volume MV <sub>spont</sub>	0 to 99 lpm	like Evita 4	0 to 99 lpm	No
Tidal volume V <sub>T</sub>	0 to 3999 ml	like Evita 4	0 to 3999 ml	Yes
spont. tidal volume V <sub>Tspont</sub>	0 to 3999 ml	like Evita 4	0 to 3999 ml	No
breathing frequency f	0 to 150 bpm	like Evita 4	0 to 80 bpm	Yes
spont. breathing frequency f <sub>spont</sub>	0 to 150 bpm	like Evita 4	0 to 80 bpm	No
breathing gas temperatur	18 to 51°C	like Evita 4	18 to 51°C	No
endexp. CO <sub>2</sub> -concentration etCO <sub>2</sub>	0 to 13.3 Vol.-%	like Evita 4 * <sup>8</sup>	-	No

\* see explanation of differences

Device	Evita 4	Evita 2 dura	Evita	Servo 300
CO <sub>2</sub> -production	0 to 999 ml/min	like Evita 4 <sup>*8</sup>	-	No
Serial dead space V <sub>ds</sub>	0 to 999 ml	like Evita 4 <sup>*8</sup>	-	No
Dead Space Ventilation V <sub>ds</sub> /VT	0 to 999 ml	like Evita 4 <sup>*8</sup>	-	-
etCO <sub>2</sub>	0 to 13.3 Vol%	like Evita 4 <sup>*8</sup>	-	-
Compliance C (calculated value)	0 to 250 ml/cmH <sub>2</sub> O	like Evita 4	0 to 250 ml/cmH <sub>2</sub> O	No
Resistance R (calculated value)	0 to 99 cmH <sub>2</sub> O/l/s	like Evita 4	0 to 99 cmH <sub>2</sub> O/l/s	No
Pulse Rate	20 to 250 bpm (option <sup>*9</sup> )	like Evita 4 <sup>*9</sup>	-	No
SpO <sub>2</sub>	0 - 100 % <sup>*9</sup>	like Evita 4 <sup>*9</sup>	-	-
<b>Displayed curves</b>				
Airway pressure Paw (t)	Yes	like Evita 4	Yes	No
Flow (t)	Yes	like Evita 4	Yes	No
Exsp. CO <sub>2</sub> -concentration	Yes	like Evita 4 <sup>*8</sup>	-	No
FCO <sub>2</sub>				
Plethmogram	Yes <sup>*9</sup>	like Evita 4	-	-

\* see explanation of differences

Device	Evita 4	Evita 2 dura	Evita	Servo 300
<b>Special Functions</b>				
Intermittend PEEP	Yes	like Evita 4	Yes	No
Inspiration Hold	Yes	like Evita 4	Yes	Yes
Apnea Ventilation	Yes	like Evita 4	Yes	Yes
Nebulization	Yes	like Evita 4	Yes	Yes
Suction Help Mode	Yes	like Evita 4	Yes	No
Switch over function	Yes	like Evita 4	Yes	Yes
<b>Physical Characteristics</b>				
Dimensions WxHxD	530 x 290 x 450 mm	530 x 290 x 420 mm	530 x 270 x 420 mm	242 x 360 240 mm (pneumatics) 427 x 110 321 mm (front panel)
Weight	27 kg	like Evita 4	27 Kg	23.2 kg
Operation Temperature	10 to 40°C	like Evita 4	15 to 40°C	
Storage Temperature	-20 to 60°C	like Evita 4	-20 to 60°C	

\* see explanation of differences

Device	Evita 4	Evita 2 dura	Evita	Servo 300
<b>Power Supply</b>				
AC-Power connection range	100V to 240 V, 50/60 Hz	like Evita 4	100V to 127 V, 220 V to 240 V, 50/60 Hz	100, 120, 220 and 240 V AC, 50/60 Hz
Power consumption	approximately 125 W	like Evita 4	approximately 125 W	40 W
DC-connection	Option, range 12 or 24V DC <sup>*10</sup>	like Evita 4	No	24 V DC
Internal battery	Option <sup>*10</sup>	like Evita 4	No	Yes
Operating time with battery	> 10 min <sup>*10</sup>	like Evita 4	No	Approx. 30 min
Typ of internal batterie	lead, free of maintenance	like Evita 4	-	lead, free of maintenance
Recharge time	< 1.5 hours	like Evita 4	-	10 hours

\* see explanation of differences

## Explanation of Differences

- 1 Evita 2 dura and Evita 4 uses the name PCV+ instead of the name BIPAP for the same ventilation mode. This pressure controlled ventilation mode is combined with free spontaneous breathing during the complete breathing cycle.
- 2 Evita 2 dura and Evita 4 uses extra setting AutoFlow for all volume controlled ventilation modes (CMV, SIMV, and MMV) instead of a separate mode (Pressure regulated volume controlled of the serve 300).
- 3 Evita 2 dura and Evita 4 distinguishes between an Adult and a Pediatric Mode, which can be selected by the user after switching on.
- 4 Evita 2 dura and Evita 4 uses instead of an I:E-adjustment the adjustment of the Inspiratory Time T<sub>insp</sub>
- 5 Evita 2 dura and Evita 4 uses the Flow Trigger for the adjustment of the trigger sensitivity in all ventilation modes
- 6 Evita 2 dura uses a different operating concept:  
11 ventilation parameter, for example VT is adjusted as follows:  
By pressing a key the parameter the adjusted parameter is selected. After that the yellow LED on the key is activated.  
Turn the central rotary knob = adjust the value  
Press the central rotary knob = confirm the value  
The yellow LED is switched off.  
The value is displayed by a 7-Segment-Display.
- 7 In Evita 2 dura for adjusting of the parameters are no potentiometer used. The adjusted value is displayed by 7-Segment-displays (each 11 parameters have an own 7-segment-display).
- 8 Evita 2 dura can be equipped with an built-in monitoring for the expiratory CO<sub>2</sub> concentration. The basic version of Evita 4 is equipped with this feature.
- 9 Evita 2 dura and Evita 4 can be equipped optional with an oxygen saturation measurement
- 10 Evita 2 dura and Evita 4 can be equipped with a DC power supply

## DESCRIPTION OF THE DEVICE:

The Evita 2 dura is a time cycled, constant volume ventilator for long term use in intensive care and inner clinical transport. The device is microprocessor controlled and provides the following ventilation modes:

- CMV/assist      Controlled Mandatory Ventilation  
Controlled and assisted volume constant ventilation with the options:
  - CPV              Continuous Positive Pressure Ventilation
  - PLV              Pressure Limited Ventilation
  - AutoFlow        For automatic regulation of inspiratory flow
  - IRV              Inversed Ratio Ventilation
- SIMV            Synchronized Intermittent Mandatory Ventilation  
Procedure for weaning patients off the ventilator after they have started spontaneous breathing. With the options:
  - PLV
  - AutoFlow
- MMV            Mandatory Minute Volume Ventilation  
Spontaneous breathing with automatic adjustment of mandatory ventilation to the patients minute volume requirement. With the options:
  - PLV
  - AutoFlow
- SB              Spontaneous Breathing  
Spontaneous breathing at ambient pressure
- CPAP            Continuous Positive Airway Pressure  
Spontaneous breathing with positive airway pressure
- ASB/PSV        Assisted Spontaneous Breathing / Pressure Support Ventilation  
Pressure supported spontaneous breathing
- PCV+            Pressure Controlled Ventilation  
Pressure controlled ventilation combined with free spontaneous breathing during the complete breathing cycle, and adjustable pressure increase to CPAP level
- APRV            Airway Pressure Release Ventilation  
Spontaneous breathing on two pressure levels with long time ranges - independently adjustable
- Apnea  
Ventilation      For switching over automatically to volume controlled mandatory ventilation if breathing stops
- ILV              Independent Lung Ventilation  
Separate, differentiated, synchronized ventilation with one Evita 2 dura for each lung

## INTENDED USE OF THE DEVICE:

The Evita 2 dura is a long term ventilator for the intensive care of patients with a body weight of at least 3 Kg up to 200 Kg.

The device is microprocessor controlled and provides the following ventilation modes:

- **CMV**                      Controlled Mandatory Ventilation
- **SIMV**                     Synchronized Intermittent Mandatory Ventilation
- **MMV**                     Mandatory Minute Volume Ventilation
- **CPAP**                    Continuous Positive Airway Pressure
- **PCV+**                    Pressure Controlled Ventilation, combined with free spontaneous breathing during the complete breathing cycle and adjustable pressure support from CPAP level.
- **APRV**                    Airway Pressure Release Ventilation  
Spontaneous breathing on two pressure levels with long time ranges - independently adjustable.
- **Apnea ventilation**      For switching over automatically to volume controlled ventilation if breathing stops.
- **ILV**                      Independent Lung Ventilation  
Separate, differentiated, synchronised ventilation with one Evita 2 dura unit for each lung.

**Application range** - Intensive Care units  
- Inner clinical transport

## SUMMARY OF THE TECHNOLOGICAL CHARACTERISTICS OF THE DEVICE COMPARED TO THE PREDICATE DEVICES:

The ventilator functions of the Evita 2 dura are substantially equivalent to the Evita 4, the Evita and the Servo 300 ventilators.

The Evita 2 dura integrates the optional functions for CO<sub>2</sub> and SpO<sub>2</sub> that are presently performed by the Evita 4.



---

Frank Clanzett  
Regulatory Affairs Specialist



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20856

Mr. Harald Kneuer  
Drager, Inc.  
Critical Care Systems  
4101 Pleasant Valley Road, Suite 100  
Chantilly, Virginia 20151

AUG 26 1997

Re: K970165  
Evita 2 Dura  
Regulatory Class: II (two)  
Product Code: 73 CBK  
Dated: May 27, 1997  
Received: May 28, 1997

Dear Mr. Kneuer:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Pre-market Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Mr. Harald Kneuer

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4648. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Thomas J. Callahan, Ph.D.  
Director  
Division of Cardiovascular,  
Respiratory, and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

