

K973696  
NOV 4 1997<sup>1</sup>

Summary of Safety & Effectiveness  
Beckman SYNCHRON® Systems Hemoglobin A1c Reagent

1.0 **Submitted By:**

Lucinda Stockert  
Senior Regulatory Specialist, Product Submissions  
Beckman Instruments, Inc.  
200 S. Kraemer Blvd., W-337  
Brea, California 92822-8000  
Telephone: (714) 961-3777  
FAX: (714) 961-4457

2.0 **Date Submitted:**

24 September 1997

3.0 **Device Name(s):**

3.1 **Proprietary Names**

SYNCHRON® Systems Hemoglobin A1c (HbA1c) Reagent

3.2 **Classification Name**

Glycosylated hemoglobin assay (21 CFR § 864.7470)

4.0 **Predicate Device(s):**

BECKMAN Reagent	Predicate	Predicate Company	Docket Number
SYNCHRON Systems Hemoglobin A1c	Tina-quant® Hemoglobin A1c Assay	Boehringer Mannheim Corporation	K934070

5.0 **Description:**

The SYNCHRON Hemoglobin A1c (HbA1c) Reagent kit contains two reagents, Hb and A1c for determining the Hemoglobin A1c concentrations in whole blood on SYNCHRON Clinical Systems.

6.0 **Intended Use:**

The SYNCHRON Hemoglobin A1c (HbA1c) Reagent kit, in conjunction with the SYNCHRON Hemoglobin A1c Calibrators, is intended for the quantitative determination of hemoglobin A1c concentration as a percentage of total hemoglobin in whole blood on SYNCHRON Systems.

7.0 **Comparison to Predicate(s):**

The following tables show similarities and differences between the predicate identified in Section 4.0 of this summary.

Reagent	Aspect/Characteristic	Comments
<b>SIMILARITIES</b>		
SYNCHRON Systems Hemoglobin A1c (HbA1c) Reagent	Intended use	Same as the predicate
	Formulation	Same source, processing as predicate
	Chemical Reaction	Same principle as the predicate
	Sample Preparation	Same as the predicate
	Calibration	Same as the predicate
	Reagent Kit Configuration	Same as the predicate
<b>DIFFERENCES</b>		
SYNCHRON Systems Hemoglobin A1c (HbA1c) Reagent	Assay parameters	SYNCHRON parameters reside in the SYNCHRON database, BMC parameters must be operator programmed
	Cartridges	Cartridges bar-coded specifically for SYNCHRON HbA1c vs non bar-coded cartridges with BMC reagent

**8.0 Summary of Performance Data:**

The data in the Premarket Notification on safety and effectiveness supports a finding of substantial equivalence to chemistry test systems already in commercial distribution. Equivalence is demonstrated through method comparison, stability, and imprecision experiments that relate results obtained from the SYNCHRON HbA1c Reagent to the Boehringer Mannheim Tina-quant® Hemoglobin A1c assay.

**Method Comparison Study Results**

Analyte	Slope	Intercept	r	Predicate Method
SYNCHRON Hemoglobin A1c (HbA1c) Reagent	0.9906	-0.19%	0.9698	Boehringer Mannheim Hemoglobin A1c Reagent on SYNCHRON CX4CE System
	0.9937	-0.25%	0.9678	Biorad Diamat™ HPLC Hemoglobin A1c Method

**Stability Study Results**

Reagent	Product Claim
SYNCHRON Hemoglobin Reagent	30 day on-board stability
SYNCHRON A1c Reagent	14 day on-board stability

**Estimated SYNCHRON HbA1c Reagent Imprecision**

Sample	Mean (%HbA1c)	S.D. (%)	%C.V.	N
<b>Within-Run Imprecision</b>				
Level 1	6.05	0.27	4.5	80
Level 2	8.75	0.26	3.0	80
Level 3	10.73	0.31	2.9	80
<b>Total Imprecision</b>				
Level 1	6.05	0.44	7.2	80
Level 2	8.75	0.52	6.0	80
Level 3	10.73	0.74	6.9	80

The Summary of Safety and Effectiveness information for the SYNCHRON Systems HbA1c Reagents are found in TAB 1 of this notice and are being submitted in accordance with the requirements of the Safe Medical Devices Act of 1990 and implementing regulation 21 CFR 807.92.



Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

NOV 24 1997

Lucinda Stockert  
Senior Regulatory Specialist  
Beckman Instruments, Inc.  
200 S. Kraemer Boulevard, M/S W-337  
P.O. Box 800  
Brea, California 92822-8000

Re: K973696  
Synchron® Systems Hemoglobin A1c (HbA1c) Reagent  
Regulatory Class: II  
Product Code: LCP  
Dated: September 24, 1997  
Received: September 26, 1997

Dear Ms. Stockert:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Pre-market Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.


Page 2

Under the Clinical Laboratory Improvement Amendments of 1988 (CLIA-88), this device may require a CLIA complexity categorization. To determine if it does, you should contact the Centers for Disease Control and Prevention (CDC) at (770) 488-7655.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Steven I. Gutman, M.D., M.B.A.  
Director  
Division of Clinical  
Laboratory Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

K973696

page 1 of 1

510(k) Number (if known): ~~Not yet assigned~~

Device Name: **SYNCHRON® Systems Hemoglobin A1c Reagent**

Indications for Use:

The SYNCHRON Systems Hemoglobin A1c (HbA1c) Reagent Kit, in conjunction with SYNCHRON HbA1c Calibrators and SYNCHRON HbA1c Hemolyzing Reagent, is intended for the quantitative determination of Hemoglobin A1c concentrations in whole blood on SYNCHRON Clinical Systems.

**21 CFR 864.7470 Glycosylated hemoglobin assay**

(a) *Identification.* A glycosylated hemoglobin assay is a device used to measure the glycosylated hemoglobins (A1a, A1b, A1c) in a patient's blood by a column chromatographic procedure. Measurement of glycosylated hemoglobin is used to assess the level of control of a patient's diabetes and to determine the proper insulin dosage for a patient. Elevated levels of glycosylated hemoglobin indicate uncontrolled diabetes in a patient.

(b) *Classification.* Class II.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

*Concurrence of CDRH, Office of Device Evaluation (ODE)*



(Division Sign-Off)

Division of Clinical Laboratory Devices

510(k) Number K973696

Prescription Use    
 (per 21 CFR 801.109)

OR

Over-the-Counter Use    
 Optional Format 1-2-96