

K981165

APR 16 1998

# "Summary of Safety & Effectiveness"

AimStep™ Combo Pregnancy is intended for the rapid qualitative identification of hCG (human Chorionic Gonadotropin) in serum or urine to aid in the early determination of pregnancy. The test utilizes a combination of monoclonal and polyclonal antibodies to selectively detect elevated levels of hCG. The assay is conducted by adding specimen to the test device and observing for the formation of colored lines. The specimen migrates via capillary action along the membrane to react with the colored conjugate. Positive specimens react with the specific antibody-hCG-colored conjugate and form a colored line in the Specimen Area of the membrane. Absence of this colored line suggests a negative result. To serve as a procedural control, a colored line in the Control Area will always appear regardless of the presence or absence of hCG.

A multi-center clinical evaluation was conducted comparing the results obtained using AimStep™ Combo Pregnancy and another commercially available serum/urine membrane test. The study included 159 urine and 73 serum specimens tested with both assays. The following results were found:


	Positive Urine Results	Negative Urine Results	
AimStep™ Combo	71	88	
Commercially Available Test	71	88	

	Positive Serum Results	Negative Serum Results	Inconclusive Serum Results
AimStep™ Combo	21	51	1
Commercially Available Test	21	51	1

AimStep™ Combo Pregnancy showed a 100% concordance with the other commercially available test.

AimStep™ Combo Pregnancy detects hCG concentrations of 20 mIU/ml and greater. The test has been standardized to the World Health Organization Third International Standard. The addition of hLH (500 mIU/ml), hFSH (1000 mIU/ml), and hTSH (1000 µIU/ml) to negative and positive serum and urine specimens showed no cross-reactivity.

  
 \_\_\_\_\_  
 Martin O'Connor, Regulatory Affairs  
 Germaine Laboratories, Inc.  
 4203 Gardendale Center, Suite 230  
 San Antonio, TX 78229 USA

3-30-98  
 \_\_\_\_\_  
 Date

Premarket Notification 510(k) Number



APR 16 1998

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Martin O'Connor  
General Manager  
Germaine Laboratories, Inc.  
4203 Gardenale Center, Suite 230  
San Antonio, Texas 78229

Re: K981165  
AimStep™ Combo Pregnancy  
Regulatory Class: II  
Product Code: JHI  
Dated: March 30, 1998  
Received: March 31, 1998

Dear Mr. O'Connor:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Pre-market Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your pre-market notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.


Page 2

Under the Clinical Laboratory Improvement Amendments of 1988 (CLIA-88), this device may require a CLIA complexity categorization. To determine if it does, you should contact the Centers for Disease Control and Prevention (CDC) at (770) 488-7655.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Steven I. Gutman, M.D., M.B.A.  
Director  
Division of Clinical  
Laboratory Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

# Indications For Use

510(k) Number: \_\_\_\_\_

Device Name: AimStep™ Combo Pregnancy

"Indications For Use" - AimStep™ Combo Pregnancy is intended for the rapid qualitative identification of hCG (human Chorionic Gonadotropin) in serum or urine to aid in the early determination of pregnancy. This test is for professional use.

(Please do not write below this point)

\_\_\_\_\_  
Concurrence of CDRH, Office of Device Evaluation (ODE)

*Loe B. Montgomery*  
\_\_\_\_\_  
(Division Sign-Off)

Division of Clinical Laboratory Devices

510(k) Number

K981165

Prescription Use              
(per 21 CFR 801.109)

or

Over-The-Counter Use