

3/12/99

K982739

<b>Summary of Safety and Effectiveness Information</b> Premarket Notification, Section 510(k)	<i>Gynex Electrodes</i> GYNEX Corporation
--	--

Regulatory Authority: Safe Medical Devices Act of 1990, 21 CFR 807.92

**1. Device Name:**

Trade Name: **Gynex Electrodes**  
 Common Name: cutting & coagulation electrodes  
 Classification Name: electrocautery, Gynecologic and accessories

**2. Establishment Name & Registration Number:**

Name: GYNEX Corporation  
 Number: 3032109

**3. Classification:**

§ 878.4120 **Electrocautery, Gynecologic and accessories.** (a) Identification. A gynecologic electrocautery is a device designed to destroy tissue with high temperatures by tissue contact with an electrically heated probe. It is used to excise cervical lesions, perform biopsies, or treat chronic cervicitis under direct visual observation. This generic type of device may include the following accessories: an electrical generator, a probe, and electrical cables. (b) Classification. Class II (performance standards).

Device Class: Class II

Classification Panel: OB/GYN

Product Code: HGI

**4. Preamendments Device:**

1. Cooper Surgical; Single-Use Tungsten Wire LEEP Electrodes, Sterile, K952483

**Gynex Corporation Device**

**Equivalent Cooper Surgical Device**

Curved electrode: 12cm x 1cm x 0.7cm

Item #R1007 ZL: 12cm x 1cm x 0.7cm

Curved electrode: 12cm x 1cm x 1cm

Item #R1010 ZL: 12cm x 1cm x 1cm

Curved electrode: 12cm x 1.5cm x 0.5cm

Item #R1505 ZL: 12cm x 1.5cm x 0.5cm

Curved electrode: 12cm x 1.5cm x 0.7cm

Item #R1507 ZL: 12cm x 1.5cm x 0.7cm

Curved electrode: 12cm x 2cm x 0.8cm

Item #R2008 ZL: 12cm x 2cm x 0.8cm

Curved electrode: 12cm x 2cm x 1cm

Item #R2010 ZL: 12cm x 2cm x 1cm

Square electrode: 12cm x 1cm x 1cm

Item #S1010 ZL: 12cm x 1cm x 1cm

Square electrode: 5.5cm x 1cm x .64cm

Item #S1007 ZL: 5.5cm x 1cm x 7.0mm

Ball electrode: 12cm x 5mm

Item #B0512 ZL: 12cm x 5mm

Ball electrode: 5.5cm x 5mm

Item #B0555 ZL: 5.5cm x 5mm

**5. Device Description:**

The *Gynex Electrodes* are made of thin 0.18mm hard tungsten wire with a high melting temperature and excellent conductiveness of electric current. The electrodes are available in three different styles intended for use in two different types of electrosurgical procedure. The three styles are curved wire loop electrodes, square wire loop electrodes and ball-type

electrodes. The two intended uses are excisional procedures (wire loop) and fulguration or coagulation procedures (ball electrode).

All electrode styles consist of an insulated conductive metallic shaft of varying length terminating in a bare metallic "ball" a "square wire loop" or a "curved wire loop". The wire loop electrodes vary in both width and depth of the loop. Electrode shaft diameter is constant at 2.5mm.

**6. Applicant Name & Address:**

GYNEX Corporation  
16700 NE 79th St., Suite 204  
Redmond, WA 98052

**7. Company Contact:**

Mr. Steve Angelo, President  
GYNEX Corporation  
16700 NE 79th St., Suite 204  
Redmond, WA 98052  
425.882.1179 \* 425.895.0115 - fax

**8. Submission Correspondent:**

Mr. David W. Schlerf  
Buckman Company, Inc.  
200 Gregory Lane, Suite C-100  
Pleasant Hill, CA 94523-3389  
925.356.2640 - 925.356.2654 - fax

**9. Voluntary Standards:**

United States Food and Drug Administration mandated performance standards for this device do not exist. Various voluntary performance standards are utilized. Voluntary standards utilized include ASTM, GYNEX Corporation. Standard Operating Procedures (SOP), vendor certification and qualification procedures, Quality Systems Regulations, ISO materials standards and ISO 9000 series quality regulations.

**10. Storage, Packaging & Sterilization Information:**

*Gynex* brand *Electrodes* are supplied "STERILE". The sterilization method is EtO gas sterilization and the exposure cycle has been validated to achieve a sterility assurance level of  $10^{-6}$ . The sterilization process, the validation procedure and the nature of the gas residuals are conducted and characterized in accordance with ANSI/AAMI/ISO 11135 -1994; EN550: 1994; AAMI method designation EOR; ISO/DIS 10993-7.2; FDA Proposed "Maximum Residue Limits and Maximum Exposure" and the Pharmacopeia of the United States.

Packaging materials are typical medical grade peel-type pouches of the generic mylar/non-woven sandwich variety. All packaging should be inspected on arrival for evidence of shipping damage. Damaged packaging may indicate the presence of unsafe product and it should not be used until carefully inspected. Any shipping damaged product should be returned promptly. Sterile product must be handled, opened and placed into use following accepted surgical sterile technique.

**11. Comparison Table:**

<b>FEATURE</b>	<b>Gynex Disposable Electrodes</b>	<b>Cooper Surgical Single-Use Wire Electrodes</b>	<b>SE?</b>
<b>Intended Use:</b>	Soft tissue cutting & coagulation Intraepithelial neoplasia Condylomata acuminata Misc. lesions; polyps, molluscum contagiosum, nevi, seborrheic keratosis, achricordon	Same	YES
<b>Design:</b>	Tugnsten wire Curved loop Square loop Brass ball electrode	Tugnsten wire Curved loop Square loop Brass ball electrode	YES
<b>Sterile:</b>	Yes	Yes	YES
<b>Sizes:</b>	<b>Curved electrode:</b> 12cm x 1cm x 0.7cm <b>Curved electrode:</b> 12cm x 1cm x 1cm <b>Curved electrode:</b> 12cm x 1.5cm x 0.5cm <b>Curved electrode:</b> 12cm x 1.5cm x 0.7cm <b>Curved electrode:</b> 12cm x 2cm x 0.8cm <b>Curved electrode:</b> 12cm x 2cm x 1cm <b>Square electrode:</b> 12cm x 1cm x 1cm <b>Square electrode:</b> 5.5cm x 1cm x 6.4mm <b>Ball electrode:</b> 12cm x 5mm <b>Ball electrode:</b> 5.5cm x 5mm  Shaft diameter 2.5mm	<b>Curved electrode:</b> 12cm x 1cm x 0.7cm <b>Curved electrode:</b> 12cm x 1cm x 1cm <b>Curved electrode:</b> 12cm x 1.5cm x 0.5cm <b>Curved electrode:</b> 12cm x 1.5cm x 0.7cm <b>Curved electrode:</b> 12cm x 2cm x 0.8cm <b>Curved electrode:</b> 12cm x 2cm x 1cm <b>Square electrode:</b> 12cm x 1cm x 1cm <b>Square electrode:</b> 5.5cm x 1cm x .7cm <b>Ball electrode:</b> 12cm x 5mm <b>Ball electrode:</b> 5.5cm x 5mm plus additional size square, ball and needle- type electrodes Shaft diameter 2.5mm	YES
<b>Material:</b>	Tungsten wire Stainless steel Brass Teflon	Tungsten wire Stainless steel Brass Teflon	YES
<b>Origin:</b>	Europe	Europe	YES
<b>Manufacturer:</b>	Gynex Corporation	Cooper Surgical	YES
<b>Product Code:</b>	HGI	HGI	YES
<b>K - Number:</b>	Pending	K952483	YES



MAR 12 1999

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850Gynex Corporation  
c/o Mr. David W. Schlerf  
Official Correspondent  
BUCKMAN Company, Inc.  
200 Gregory Lane, Suite C-100  
Pleasant Hill, CA 94523-3389Re: K982739  
GYNEX Electrodes  
Dated: November 20, 1998  
Received: February 24, 1999  
Regulatory Class: II  
21 CFR 878.4120/Procode: 85 HGI

Dear Mr. Schlerf:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4613. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597, or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

CAPT Daniel G. Schultz, M.D.  
Acting Director, Division of Reproductive,  
Abdominal, Ear, Nose and Throat,  
and Radiological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number : K982739

Device Name(s):

**GYNEX Electrodes**


**Intended Use of the Device:**

1. Soft tissue cutting & coagulation
2. Intraepithelial neoplasia
3. Condylomata acuminata
4. Misc. lesions; polyps, molluscum contagiosum, nevi, seborrheic keratosis, achricordon, etc.

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NECESSARY

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
David A. Segerson  
(Division Sign-Off)  
Division of Reproductive, Abdominal, ENT,  
and Radiological Devices  
510(k) Number K982739 / S<sup>001</sup>

Prescription Use X  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use \_\_\_\_\_  
(Optional format 1-2-96)