

JAN 22 1999

K984321

Chapter 1 – Summary Information

510(k) Summary

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92

The assigned 510(k) number is: _____.

1. Submitter name, address, contact

Ortho-Clinical Diagnostics, Inc.
100 Indigo Creek Drive
Rochester, New York 14626-5101
(716) 453-3607

Contact Person: Anne Zavertrnik

Date 510(k) prepared: November 30th, 1998

2. Device Name

Vitamin B12 assay

Trade or Proprietary Name: VITROS Immunodiagnostic Products Vitamin B12 assay

Common Name: Vitamin B12 assay

Classification Name: Vitamin B12 assay for the *in vitro* quantitative measurement of vitamin B12 in human serum and plasma (EDTA or heparin).

3. Predicate Device

The VITROS Immunodiagnostic Products Vitamin B12 assay is substantially equivalent to Bio-Rad Quantaphase II B₁₂ Radioassay.

510(k) Summary, continued

4. Device Description

The VITROS Immunodiagnostic System uses luminescence as the signal in the quantitative and semi-quantitative determination of selected analytes in human body fluids, commonly serum, plasma and urine. Coated microwells are used as the solid phase separation system.

The system is comprised of three main elements:

1. The VITROS Immunodiagnostic Products range of products, in this case VITROS Immunodiagnostic Products Vitamin B12 Reagent Pack 1/2, VITROS Immunodiagnostic Products Vitamin B12/Folate Reagent Pack 3, VITROS Immunodiagnostic Products Vitamin B12 Calibrators and the VITROS Immunodiagnostic System.
2. The VITROS Immunodiagnostic System - instrumentation, which provides automated use of the immunoassay kits. The VITROS Immunodiagnostic System was cleared for market by a separate 510(k) pre-market notification (K962919).
3. Common reagents used by the VITROS System in each assay. The VITROS Immunodiagnostic Products Signal Reagent and VITROS Immunodiagnostic Products Universal Wash Reagent were cleared as part of the VITROS Immunodiagnostic Products Total T3 510(k) pre-market notification (K984310).

The VITROS System and common reagents are dedicated specifically only for use with the VITROS Immunodiagnostic Products range of immunoassay products.

5. Device Intended Use

The VITROS Vitamin B12 assay is intended for the *in vitro* quantitative measurement of vitamin B12 in human serum and plasma (EDTA or heparin), to aid in the differential diagnosis of anemia.

6. Comparison to Predicate Device

The VITROS Immunodiagnostic Products Vitamin B12 assay is substantially equivalent to Bio-Rad Quantaphase II B₁₂ Radioassay (predicate device), which was cleared by FDA (K935286) for IVD use.

The relationship between the VITROS Vitamin B12 assay and the predicate device, determined by Deming's Regression, is:

VITROS Vitamin B12 assay = $0.984 \times \text{Bio-Rad Quantaphase II B}_{12} \text{ Radioassay} + 9.59$ (pg/mL).

Comparisons of the VITROS Vitamin B12 assay and the predicate device were performed with samples from a variety of clinical categories.

510(k) Summary, continued

In addition to the studies mentioned above, tests were performed to obtain analytical sensitivity, specificity, precision, dilution and expected values. Refer to the VITROS Vitamin B12 assay package insert for VITROS Vitamin B12 assay results.

Table 1 lists the similarities and differences of the device characteristics between the VITROS Vitamin B12 assay with the predicate device, Bio-Rad Quantaphase II B₁₂ Radioassay.

Table 1 List of the assay characteristics

Device Characteristic	VITROS Vitamin B12 assay	Predicate Device
Calibration range	0 – 2000 pg/mL	0 – 2000 pg/mL
Basic principle	Solid phase immunoassay	Radioassay
Tracer	Enzyme labeled	⁵⁷ Co
Instrumentation	VITROS Immunodiagnostic System	Gamma Counter
Sample type	Serum, plasma (EDTA or heparin)	Serum, plasma (EDTA)
Sample volume	30 µL	200 µL
Incubation time and temperature	58 minutes at 37° C	1 hour at room temperature

7. Conclusions

The data presented in the pre-market notification demonstrate that the VITROS Vitamin B12 assay performs substantially equivalent to the predicate device, which was cleared by FDA (K935286) for IVD use.

Equivalence was demonstrated using currently commercially available reagents along with patient specimens covering a variety of clinical categories.

The data presented in the premarket notification provide a reasonable assurance that the VITROS Vitamin B12 assay is safe and effective for the stated intended use.


Page 2

Under the Clinical Laboratory Improvement Amendments of 1988 (CLIA-88), this device may require a CLIA complexity categorization. To determine if it does, you should contact the Centers for Disease Control and Prevention (CDC) at (770) 488-7655.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification"(21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597, or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink that reads "Steven Gutman". The signature is written in a cursive style.

Steven I. Gutman, M.D, M.B.A.
Director
Division of Clinical
Laboratory Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Statement of Intended Use

510(k) Number (if known):

K984321

Device Name:

1. VITROS Immunodiagnostic Products Vitamin B12 Reagent Pack 1/2
2. VITROS Immunodiagnostic Products Vitamin B12/Folate Reagent Pack 3
3. VITROS Immunodiagnostic Products Vitamin B12 Calibrators

Indications for Use:

1 & 2. The VITROS Immunodiagnostic Products Vitamin B12 Reagent Pack 1/2 and the VITROS Immunodiagnostic Products Vitamin B12/Folate Reagent Pack 3 – for the *in vitro* quantitative measurement of vitamin B12 in human serum and plasma (EDTA or heparin) to aid in the differential diagnosis of anemia.

3. The VITROS Immunodiagnostic Products Vitamin B12 Calibrators – for *in vitro* use in the calibration of the VITROS Immunodiagnostic System for the quantitative measurement of vitamin B12 in human serum and plasma (EDTA or heparin).

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓
(Per 21 CFR 801.109)

OR

Over-The-Counter Use _____

(Optional Format 1-2-96)

Carol Benam for Jean Cooper
(Division Sign-Off)

Division of Clinical Laboratory Devices

510(k) Number K984321