



March 17, 1999

K990927

11-12-99

510(k) SUMMARY
Page 1 of 3

**Food and Drug Administration
Center for Devices and Radiological Health
Document Mail Center (HFZ-401)
9200 Corporate Blvd.
Rockville, Maryland 20850**

Attention: Document Control Clerk

Subject: **SIGMA 3400 Monoplace Hyperbaric Chamber Summary of Safety & Effectiveness Information**

Gentlemen:

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA of 1990 and 21 CFR 807.92.

Substantial Equivalence: The Perry Baromedical Corporation **SIGMA 3400** Monoplace Hyperbaric Chamber is substantially equivalent to the Perry Baromedical Corporation **SIGMA I** Monoplace Hyperbaric Chamber, [510(k) K832127], and the Perry Baromedical Corporation **SIGMA Plus** Monoplace Hyperbaric Chamber, [510(k) 974863].

The **SIGMA 3400** has the same intended use and the same technological characteristics as the **SIGMA I** and **SIGMA Plus** predicate devices. The scope and operation of the **SIGMA 3400** is identical to that of the **SIGMA I** and **SIGMA Plus** in the following areas:

- 1 The critical component of both systems consists of an ASME pressure vessel that is designed, fabricated, and tested in accordance with the requirements of the ASME Boiler and Pressure Vessel Code, Section VIII, Division 1; and ASME PVHO-1, Safety Standard for Pressure Vessels for Human Occupancy.
- 2 The pressure vessel consists of certified metal components constructed of materials in conformance with the ASME Boiler and Pressure Vessel Code, Section II - Materials. In addition to the metal parts, the pressure retaining boundary consists of an acrylic plastic cylinder, designed and constructed to the requirements of ASME PVHO-1, Safety Standard for Pressure Vessels for Human Occupancy.

- 3 The pressure vessel is protected from accidental over pressurization by an ASME "UV" stamped pressure relief valve.
- 4 The medium used for pressurization of the chamber is medical grade oxygen.
- 5 Treatment capacity of the chamber is one patient.
- 6 Communication between the patient / attendants is through an intrinsically safe communications system.
- 7 The **SIGMA I**, **SIGMA Plus** and **SIGMA 3400** chambers are used to provide Hyperbaric Oxygen Therapy, prescribed by a licensed physician, where medical grade oxygen is administered to patients while under pressure, by the patient breathing oxygen from the chamber atmosphere that is pressurized with pure oxygen. The breathing of oxygen at elevated pressure promotes the movement of oxygen into the patient's tissues.
- 8 The chamber pressurization and depressurization control system is completely pneumatic, with no electrical or electronic components.
- 9 Chambers are Class "B" monoplace hyperbaric systems.

Differences do exist between the **SIGMA 3400**, the **SIGMA Plus** and the **SIGMA I** chambers; they include:

- 1 The PVHO-1 acrylic cylinder in the **SIGMA 3400** chamber is 33.3" inside diameter and 90" long; the PVHO-1 acrylic cylinder in the **SIGMA Plus** chamber is 40" inside diameter by 65" long, and the PVHO-1 acrylic cylinder in the **SIGMA I** chamber is 25.25" inside diameter by 80" long.
2. The acrylic cylinder of the **SIGMA 3400** has been designed for a operational temperature range of 32 degrees F to 100 degrees F. The **SIGMA I** and **SIGMA Plus** cylinders were designed for a 32 degree F to 125 degree F operating range. Based on the actual operating environment of a hospital, the higher 125 degree temperature limit would never be reached in normal operation, as the chamber rooms are required to be maintained at temperatures less than 75 degrees F for patient comfort.
3. The pneumatic controls for the **SIGMA 3400** chamber and the **SIGMA I** chamber are mounted in a panel mounted on one side of the chamber. The pneumatic controls for the **SIGMA Plus** chamber are mounted in a control box that is suspended from a moveable arm.


- 4 The **SIGMA 3400** monoplace chamber has an aluminum chassis with integral side covers; the **SIGMA Plus** Chamber has gelcoated fiberglass cosmetic covers on the sides and ends of the chamber; and the **SIGMA I** Hyperbaric Chamber has aluminum side cosmetic covers as part of the chassis.

5. The **SIGMA 3400** monoplace chamber has provision for changing from oxygen to air pressurization during treatment, which allows the operator more control of the environment inside the chamber. The **SIGMA Plus** and **SIGMA I** chambers are oxygen only pressurization.

The changes to the **SIGMA I** and **SIGMA Plus** systems resulting in the **SIGMA 3400** Hyperbaric Chamber enhance its ease of operation, and do not adversely effect the safety and effectiveness of the device. **The basic treatment method is identical, only the chamber acrylic window diameter, length and operating temperature range has been changed**

None of the above information is confidential and all may be made available to the public upon written request.

Perry Baromedical Corporation


James P. Dodson
Quality Assurance Manager

3-17-99
Date



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

NOV 12 1999

Mr. James P. Dodson
Perry Baromedical Corp.
3660 Interstate Parkway
Riviera Beach, FL 33404-3411

Re: K990927
SIGMA 3400 Monoplace Hyperbaric Chamber
Regulatory Class: II (two)
Product Code: 73 CBF
Dated: August 27, 1999
Received: August 30, 1999

Dear Mr. Dodson:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Pre-market Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Mr. James P. Dodson

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4648. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597, or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Celia M. Witten, Ph.D., M.D.
Acting Director
Division of Cardiovascular,
Respiratory, and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure



**PERRY
BAROMEDICAL**

3660 Interstate Parkway
Riviera Beach
Florida 33404 USA
800•741•4376
561•840•0395
FAX: 561•840•0398

510(k) K 990927

SIGMA 3400 Monoplace Hyperbaric Chamber

INDICATIONS FOR USE

The **SIGMA 3400 Monoplace Hyperbaric Chamber** is indicated for use for the following clinical medical conditions in accordance with guidelines established by the Undersea and Hyperbaric Medical Society, as follows:


1. Air or Gas Embolism
2. Carbon Monoxide / Smoke Inhalation
3. Compromised Skin Grafts and Flaps
4. Crush Injuries / Acute Traumatic Ischemias
5. Decompression Sickness
6. Enhanced Healing in Selected Problem Wounds
7. Exceptional Blood Loss (Anemia)
8. Gas Gangrene (Clostridial Myonecrosis)
9. Intracranial Abscess
10. Necrotizing Soft Tissue Infections
11. Radiation Tissue Damage (Osteoradionecrosis)
12. Refractory Osteomyelitis
13. Thermal Burns

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use
(Per 21 CFR 801.109)

OR

Over-The-Counter-Use
(Optional Format 1-2-96)


(Division Sign-Off)
Division of Cardiovascular, Respiratory,
and Neurological Devices
510(k) Number K990927