

JUN 9 1999

K991255

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

CADD-TPN[®] Administration Set with 1.2 μ Filter and Add-on Anti-siphon Valve

April 9, 1999

I. GENERAL INFORMATION

Applicant's Name and Address: SIMS Deltec, Inc.
1265 Grey Fox Road
St. Paul, MN 55112

Contact Person: Edward W. Numainville
Vice President, Regulatory Affairs and
Quality Systems
Tel. (651) 628-7166

Common/Usual Name: Administration Set

Proprietary Name: CADD-TPN[®] Administration Set with 1.2 μ
Filter and Add-on Anti-siphon Valve

Equivalence Device Comparison: CADD-TPN[®] Administration Set with 1.2 μ
Filter

CADD[®] Extension Set with Anti-siphon Valve

CADD[®] Administration Set with TOTM
Plasticized Tubing

II. DEVICE DESCRIPTION

The purpose of this submission is to offer an alternate CADD-TPN[®] Administration Set with 1.2 μ Filter and Add-on Anti-siphon Valve. The set is made up of the following components: bag spike, tubing, set plate, 1.2 μ filter, tubing clamp, injection site, male luer connector with protective cap, and anti-siphon valve. The anti-siphon valve is designed to protect against unregulated gravity infusion that can result from an improperly attached set plate to the CADD-TPN[®] Ambulatory Infusion Pump.

III. INTENDED USE OF THE DEVICE

The CADD-TPN[®] Administration Set is designed for use with the CADD-TPN[®] pump for the administration of nutritional and other solutions or fluids intravenously from an IV bag.

IV. DEVICE COMPARISON

	CADD-TPN® Administration Set with 1.2 µ Filter and Add-on Anti-siphon Valve	CADD-TPN® Administration Set with 1.2 µ Filter	CADD® Extension Set with Anti-siphon Valve	CADD® Administration Set with TOTM Plasticized Tubing
MANUFACTURER	SIMS Deltec, Inc.	SIMS Deltec, Inc.	SIMS Deltec, Inc.	SIMS Deltec, Inc.
510(K) NUMBER	Subject Device	K902374	K942046	K933390
INDICATIONS FOR USE	The CADD-TPN® Administration Set is designed for use with the CADD-TPN® pump for the administration of nutritional and other solutions or fluids intravenously from an IV bag.	The CADD-TPN® Administration Set is designed for use with the CADD-TPN® pump for the administration of nutritional and other solutions or fluids intravenously from an IV bag.	For use with the MEDICATION CASSETTE™ Reservoir.	The CADD® Administration Set is designed to be used with a variety of Deltec CADD® Pumps to allow fluid delivery from a remote bag.
DIMENSIONS (Nominal)				
LENGTH	114 in.	100 in.	30, 45, or 60 in.	77 in.
I.D.	0.100 in.	0.100 in.	0.039 in.	0.100 in.
O.D.	0.164 in.	0.164 in.	0.106 in.	0.164 in.
FILTER	YES	YES	NO	NO
ANTI-SIPHON VALVE	YES	NO	YES	YES

V. SUMMARY OF STUDIES

A. **Functional Testing**

Functional specification testing was performed on the alternate filter and anti-siphon valve.

Biocompatibility testing was conducted.

B. **Clinical Studies**

Clinical studies were not deemed necessary regarding the new set due to its similarity in materials, design and function to current commercially available CADD® Administration Sets and Extension Sets.

C. Conclusions Drawn from the Studies

The results of the testing indicated that the filter functions according to specification and the set meets the biocompatibility requirements. Therefore, the CADD-TPN[®] Administration Set with 1.2 μ Filter and Add-on Anti-siphon Valve is considered acceptable for human use.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville, MD 20850

JUN 9 1999

Mr. Edward W. Numainville
Vice President
Regulatory Affairs and Quality Systems
SIMS Deltec, Incorporated
1265 Grey Fox Road
St. Paul, Minnesota 55112

Re: K991255
Trade Name: CADD-TPN® Administration Set with 1.2 μ
Filter and Add-on Anti-siphon Valve
Regulatory Class: II
Product Code: FPA
Dated: April 9, 1999
Received: April 13, 1999

Dear Mr. Numainville:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Pre-market Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Page 2 - Mr. Numainville

Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski
Director
Division of Dental, Infection Control,
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

K 991255

510(k) Number (if known): _____

Device Name: CADD-TPN[®] Administration Set with 1.2 μ Filter and Add-on Anti-siphon Valve

Indications for Use:

“The CADD-TPN[®] Administration Set is designed for use with the CADD-TPN[®] pump for the administration of nutritional and other solutions or fluids intravenously from an IV bag.”

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

OR

Over-The Counter Use _____

Peterson Cuervo
(Division Sign-Off)
Division of **Dental, Infection Control,**
and General **Hospital Devices**
510(k) Number K 991255