

FDA/CDRH Premarket 510(k) Notification  
 GE Medical Systems - LOGIQ 200 MD  
 May 24, 1999

## Section 2:

*Summary of Safety and Effectiveness  
 Prepared in accordance with 21 CFR Part 807.92(c).*



GE Medical Systems

General Electric Company  
 P.O. Box 414, Milwaukee, WI 53201

### Section a):

1. Submitter: GE Medical Systems  
 PO Box 414  
 Milwaukee, WI 53201

Contact Person: Allen Schuh,  
 Manager, Safety and Regulatory Engineering  
 Telephone: 414-647-4385, Fax: 414-647-4090

Date Prepared: May 24, 1999

2. Device Name: GE LOGIQ 200 MD Diagnostic Ultrasound System  
 Ultrasonic Pulsed Echo Imaging System, 21 CFR 892.1560, 90-IYO

3. Marketed Device: GE LOGIQ  $\alpha$ 200 Diagnostic Ultrasound System - K960700, 28 AUG 1996

4. Device Description: The GE LOGIQ 200 MD consists of a mobile console approximately 40 cm wide, 60 cm deep and 125 cm high, weighing approximately 75 kg. The user interface consists of a keyboard and monochrome video monitor. The system is designed for use in B, M and B/M scanning modes and supports linear, convex, micro convex and phased array probes.

5. Indications for Use: The GE LOGIQ 200 MD is a general purpose imaging system intended for use in obstetrics, gynecology, urology and general radiology examinations by a qualified physician to aid in the diagnosis and evaluation of soft tissues by generating two dimensional images, time-motion images and biometric data. Specific Clinical applications included in the indications for use are: Fetal, Abdominal, Intraoperative (abdominal), Pediatric, Small organs (breast, neck, chest, male and female reproductive organs, limbs & extremities), Adult cephalic, Neonatal cephalic, Cardiac (adult & pediatric), Transvaginal, Transrectal and Urological.

6. Comparison with Predicate Device: The GE LOGIQ 200 MD Diagnostic Ultrasound System is of a comparable type and substantially equivalent to the currently marketed GE LOGIQ  $\alpha$ 200. It has similar technological characteristics, is comparable in key safety and effectiveness features, uses the same general design, construction, and materials, and has the same intended uses and operating modes as the predicate device.

### Section b):

1. Non-clinical Tests: The device has been evaluated for acoustic output, biocompatibility, and thermal, electrical and mechanical safety, and has been found to conform with applicable medical device safety standards. The ultrasound generator modifications require re-verification of the acoustic output levels which will be provided in the 510(k) special report.

2. Clinical Tests: None required.

3. Conclusion: Intended uses and other key features are consistent with traditional clinical practice, FDA guidelines, and established methods of patient examination. The design and development process of the manufacturer conforms with 21 CFR 820 Quality System Regulation and ISO 9001 & 13485 quality system standards. The product is designed to conform with applicable medical device safety standards and compliance is verified through independent evaluation with ongoing production surveillance. Diagnostic ultrasound has accumulated a long history of safe and effective performance. Therefore, it is the opinion of GE Medical Systems that the GE LOGIQ 200 MD is substantially equivalent with respect to safety and effectiveness to devices currently cleared for market.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

JUL 12 1999

General Electric Medical Systems  
c/o Chantel Carson  
Engineering Services  
Underwriters Laboratories, Inc.  
353 Pfingsten Road  
Northbrook, IL 60062-2096

Re: K992208  
Trade Name: GE LOGIQ 200 MD Diagnostic Ultrasound System  
Regulatory Class: II  
Product Code: 90-IYO and 90-ITX  
Dated: June 22, 1999  
Received: June 30, 1999

Dear Ms. Carson:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the GE LOGIQ 200 MD Diagnostic Ultrasound System, as described in your premarket notification:

Transducer Model Numbers

CBF, CAE, MTZ, CZB, LH, LE, LI, LT, LB, LD, CS, 9Lb, SY, ERB, and 3Cb

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval) it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic QS inspections, the FDA will verify such assumptions. Failure to comply with the GMP regulation may result in

regulatory action. In addition, the Food and Drug Administration (FDA) may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification does not affect any obligation you may have under sections 531 and 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Please be advised that the determination above is based on the fact that no medical devices have been demonstrated to be safe and effective for in vitro fertilization or percutaneous umbilical blood sampling, nor have any devices been marketed for these uses in interstate commerce prior to May 28, 1976, or reclassified into class I (General Controls) or class II (Special Controls). FDA considers devices specifically intended for in vitro fertilization and percutaneous umbilical blood sampling to be investigational, and subject to the provision of the investigational device exemptions (IDE) regulations, 21 CFR, Part 812. Therefore, your product labeling must be consistent with FDA's position on this use.

This determination of substantial equivalence is granted on the condition that prior to shipping the first device, you submit a postclearance special report. This report should contain complete information, including acoustic output measurements based on production line devices, requested in Appendix G, (enclosed) of the Center's September 30, 1997 "Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers." If the special report is incomplete or contains unacceptable values (e.g., acoustic output greater than approved levels), then the 510(k) clearance may not apply to the production units which as a result may be considered adulterated or misbranded.

The special report should reference the manufacturer's 510(k) number. It should be clearly and prominently marked "ADD-TO-FILE" and should be submitted in duplicate to:


Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center (HFZ-401)  
9200 Corporate Boulevard  
Rockville, Maryland 20850

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4591. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

If you have any questions regarding the content of this letter, please contact Robert Phillips, Ph.D. at (301) 594-1212.

Sincerely yours,



CAPT Daniel G. Schultz, M.D.  
Acting Director  
Division of Reproductive, Abdominal,  
Ear, Nose, and Throat, and  
Radiology Devices  
Office of Device Evaluation  
Center for Devices and Radiological  
Health

Enclosure(s)

Page 4 - Ms. Chantel Carson

cc:  
HFZ-401  
HFZ-470  
D.O.

draft:  
final:RPhillips:7.9.99

**Diagnostic Ultrasound Indications for Use Form**

**GE LOGIQ 200 MD System**

Intended Use: 2D ultrasound imaging or motion analysis as follows:


Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		P	P						P	
Abdominal		P	P						P	
Intraoperative (specify)		P	P						P	
Intraoperative Neurological										
Pediatric										
Small Organ (specify)		P	P						P	
Neonatal Cephalic		P	P						P	
Adult Cephalic		P	P						P	
Cardiac		P	P						P	
Transesophageal										
Transrectal		P	P						P	
Transvaginal		P	P						P	
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)		P	P						P	

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Cardiac is Adult and Pediatric. Small Organ includes: breast, neck, chest, male and female reproductive organs, limbs, & extremities. Other includes urological. Intraoperative includes abdominal organs. Combined mode is B/M. Previous market clearance: K960700

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal, ENT,  
 and Radiological Devices

510(k) Number K992208

  
 Prescription User (Per 21 CFR 801.109)

**Diagnostic Ultrasound Indications for Use Form**  
**GE LOGIQ 200 MD with CBF Transducer**

Intended Use: 2D ultrasound imaging or motion analysis:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		P	P						P	
Abdominal		P	P						P	
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Combined mode is B/M.

Previous market clearance: K960700

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)  
 Division of Reproductive, Abdominal, ENT,  
 and Radiological Devices

510(k) Number K992208

Prescription Use (Per 21 CFR 801.109)

**Diagnostic Ultrasound Indications for Use Form**  
**GE LOGIQ 200 MD with CAE Transducer**

Intended Use: 2D ultrasound imaging or motion analysis:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		P	P						P	
Abdominal		P	P						P	
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac		P	P						P	
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										


N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Combined mode is B/M. Cardiac includes adult and pediatric.

Previous market clearance: K960700

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal, ENT,  
 and Radiological Devices

510(k) Number K992208

Prescription User (Per 21 CFR 801.109)

**Diagnostic Ultrasound Indications for Use Form**  
**GE LOGIQ 200 MD with MTZ Transducer**

Intended Use: 2D ultrasound imaging or motion analysis:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		P	P						P	
Abdominal		P	P						P	
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal		P	P						P	
Transvaginal		P	P						P	
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)		P	P						P	

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Combined mode is B/M. Other includes urological.

Previous market clearance: K960700

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)  
 Division of Reproductive, Abdominal, ENT,  
 and Radiological Devices

510(k) Number K992208

Prescription User (Per 21 CFR 801.109)

**Diagnostic Ultrasound Indications for Use Form**  
**GE LOGIQ 200 MD with CZB Transducer**

Intended Use: 2D ultrasound imaging or motion analysis:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)		P	P						P	
Neonatal Cephalic		P	P						P	
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Combined mode is B/M.

Previous market clearance: K960700

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)  
 Division of Reproductive, Abdominal, ENT,  
 and Radiological Devices

510(k) Number K992208

Prescription User (Per 21 CFR 801.109)

**Diagnostic Ultrasound Indications for Use Form**  
**GE LOGIQ 200 MD with LH Transducer**

Intended Use: 2D ultrasound imaging or motion analysis:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)		P	P						P	
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										

N= new indication; P= previously cleared by FDA; E= added under Appendix E


Additional Comments: Small Organ includes: breast, neck, chest, male and female reproductive organs, limbs, & extremities. Combined mode is B/M. Previous market clearance: K960700

\_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
 Prescription User (Per 21 CFR 801.109)

  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal, ENT,  
 and Radiological Devices  
 510(k) Number K992208

**Diagnostic Ultrasound Indications for Use Form**  
**GE LOGIQ 200 MD with LE Transducer**

Intended Use: 2D ultrasound imaging or motion analysis:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		P	P						P	
Abdominal		P	P						P	
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										


N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Combined mode is B/M.

Previous market clearance: K960700

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal, ENT,  
 and Radiological Devices

510(k) Number K 99 2208

✓  
 Prescription User (Per 21 CFR 801.109)

**Diagnostic Ultrasound Indications for Use Form**  
**GE LOGIQ 200 MD with LI Transducer**

Intended Use: 2D ultrasound imaging or motion analysis:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal		P	P						P	
Intraoperative (specify)		P	P						P	
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Combined mode is B/M.

Previous market clearance: K960700

Intraoperative includes abdominal organs.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)

Division of Reproductive, Abdominal, ENT,  
 and Radiological Devices

510(k) Number K992208

✓  
 Prescription User (Per 21 CFR 801.109)

**Diagnostic Ultrasound Indications for Use Form**  
**GE LOGIQ 200 MD with LT Transducer**

Intended Use: 2D ultrasound imaging or motion analysis:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal		P	P						P	
Intraoperative (specify)		P	P						P	
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Combined mode is B/M.

Previous market clearance: K960700

Intraoperative includes abdominal organs.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)

Division of Reproductive, Abdominal, ENT,  
 and Radiological Devices

510(k) Number K992208

Prescription User (Per 21 CFR 801.109)

**Diagnostic Ultrasound Indications for Use Form**  
**GE LOGIQ 200 MD with LB Transducer**

Intended Use: 2D ultrasound imaging or motion analysis:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		P	P						P	
Abdominal		P	P						P	
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										

N= new indication; P= previously cleared by FDA; E= added under Appendix E


Additional Comments: Combined mode is B/M.

Previous market clearance: K960700

\_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)  
 Division of Reproductive, Abdominal, ENT,  
 and Radiological Devices

510(k) Number K99 2208

✓  
 Prescription User (Per 21 CFR 801.109)

**Diagnostic Ultrasound Indications for Use Form**  
**GE LOGIQ 200 MD with LD Transducer**

Intended Use: 2D ultrasound imaging or motion analysis:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal		P	P						P	
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Combined mode is B/M.

Previous market clearance: K960700

\_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

**Concurrence of CDRH, Office of Device Evaluation (ODE)**

✓  
 Prescription User (Per 21 CFR 801.109)

  
 \_\_\_\_\_  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal, ENT,  
 and Radiological Devices  
 510(k) Number K 99 22 08

**Diagnostic Ultrasound Indications for Use Form**  
**GE LOGIQ 200 MD with CS Transducer**

Intended Use: 2D ultrasound imaging or motion analysis:


Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		N	N						N	
Abdominal		N	N						N	
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac		N	N						N	
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Cardiac is Adult and Pediatric. Combined mode is B/M.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal, ENT,  
 and Radiological Devices  
 510(k) Number K992209

  
 Prescription User (Per 21 CFR 801.109)

**Diagnostic Ultrasound Indications for Use Form**  
**GE LOGIQ 200 MD with 9Lb Transducer**

Intended Use: 2D ultrasound imaging or motion analysis:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)		N	N						N	
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Combined mode is B/M. Small organ includes: breast, neck, chest, male and female reproductive organs, limbs, & extremities.

---



---



---



---

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)  
 Division of Reproductive, Abdominal, ENT,  
 and Radiological Devices

510(k) Number   K99 2208  

✓  
 Prescription User (Per 21 CFR 801.109)

**Diagnostic Ultrasound Indications for Use Form**  
**GE LOGIQ 200 MD with SY Transducer**

Intended Use: 2D ultrasound imaging or motion analysis:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		N	N						N	
Abdominal		N	N						N	
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac		N	N						N	
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										

N= new indication; P= previously cleared by FDA; E= added under Appendix E


Additional Comments: Combined mode is B/M. Small organ includes breast, testes, thyroid.

Cardiac includes adult and pediatric.

\_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)  
 Division of Reproductive, Abdominal, ENT,  
 and Radiological Devices

510(k) Number 1K992208

Prescription User (Per 21 CFR 801.109)

**Diagnostic Ultrasound Indications for Use Form**  
**GE LOGIQ 200 MD with ERB Transducer**

Intended Use: 2D ultrasound imaging or motion analysis:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal										
Abdominal		N	N						N	
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal		N	N						N	
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)		N	N						N	

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Combined mode is B/M.


Other includes: Urological

\_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

✓  
 Prescription User (Per 21 CFR 801.109)

  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal, ENT,  
 and Radiological Devices  
 510(k) Number K992208

**Diagnostic Ultrasound Indications for Use Form**  
**GE LOGIQ 200 MD with 3Cb Transducer**

Intended Use: 2D ultrasound imaging or motion analysis:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic										
Fetal		N	N						N	
Abdominal		N	N						N	
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										


N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments: Combined mode is B/M.

\_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

**Concurrence of CDRH, Office of Device Evaluation (ODE)**



**(Division Sign-Off)**  
 Division of Reproductive, Abdominal, ENT,  
 and Radiological Devices

510(k) Number 1K992208

✓  
 Prescription User (Per 21 CFR 801.109)