

JUN 21 2000

K992606

## 510(k) Summary

**Prepared:** July 15, 1999

**Submitter:**

Company Name: Canon USA, Inc. (U.S. designated agent for Canon Inc.)  
Company Address: One Canon Plaza  
Lake Success, NY 11042  
Contact Person: Ken Shadoff, Senior Product Safety Engineer  
Phone Number: (516) 328-5602  
Fax Number: (516) 328-5169

**Proposed Device:**

Reason For 510(k): New Model  
Manufacturer: Canon Inc.  
Trade Name: Canon Laser Blood Flowmeter  
Model Name: CLBF Model 100  
Classification Name: HLI, Ophthalmoscope, AC-Powered  
FDA 510(k) #: To be assigned

**Predicate Device:**

Manufacturer: Heidelberg Engineering GmbH  
Trade Name: Heidelberg  
Model Name: Retina Flowmeter  
Classification Name: HLI, Ophthalmoscope, AC-Powered  
FDA 510(k) #: K943955

**Description Of Device:**

The Canon Laser Blood Flowmeter (CLBF) Model 100 allows for fundus observation with velocity measurement of retinal blood flow.

**Intended Use:**

This device is intended for use in observing patient's fundus (retina) with the information of blood flow (velocity) in a retinal vessel. The information of blood flow is also obtained by CLBF Model 100, and it is intended to help doctor's reading of retinal images.

**Technical Characteristics:**

Please refer to the attached COMPARISON CHART.

**Non-clinical Tests:**

Software evaluation has been performed to support a claim of substantial equivalence. Documented non-clinical studies have been included to support a claim of substantial equivalence.

**Clinical Tests:**

A published clinical study was included to support a claim of substantial equivalence.

Comparison with the similar product of the other company (1/2)

Approved product : Heidelberg Engineering, Heidelberg Retina Flowmeter

July 15, 1999

Manufacturer		Canon Inc.		Heidelberg Engineering	
product name		Canon Laser Blood Flowmeter CLBF model 100		Heidelberg Retina Flowmeter	
Intended use		Fundus observation with retinal blood flow measurement			
specification / function	fundus observation	method	similar to conventional fundus camera with eye drop ( through optical finder / on PC monitor (optional ) )	laser scanning ophthalmoscope with no eye drop ( displayed on PC monitor )	
		observation angle	30		max. 20x 20
		illumination light	560 - 650 nm Halogen lamp		670 nm / 100 uW laser diode
		diopter compensation range for examinee's eye	-10 - +10 Diopter or more		-12 - +12 Diopter
	velocity measurement	principle	laser Doppler		laser Doppler
		measurement laser	675 nm / max. 300uW ( usually 200uW ) laser diode		785 nm / 100uW laser diode
		target	blood flow in retinal vessels		blood flow in retinal vessels

Comparison with the similar product of the other company (2/2)

Product Name		Canon Laser Blood Flowmeter CLBF model 100	Heidelberg Retina Flowmeter
specification / function	fundus observation	2 sec for data acquisition	2 sec
	measurement duration	given in a real dimension i.e. velocity : [mm/sec], flow rate : [ul/min]	given as relative value
	measurement result	numerically : velocity flow rate vessel diameter graphically : time variation of velocity, flow rate	numerically : velocity, flow, blood volume graphically : perfusion map of retina velocity / flow / volume
	others	retina (one or both eyes) general population limit exposure to laser light UL 544, IEC 825 pending 510(k)	same same same unknown
	anatomical regions	retina (one or both eyes)	same
	target population	general population	same
	physical safety	limit exposure to laser light UL 544, IEC 825	same
	compliance with standards	UL 544, IEC 825 pending 510(k)	unknown
	status of device	pending 510(k)	post-amendment device, K943955
	dimensions ( W x D x H )	main unit : 320 x 560 x 620 power control unit : 230 x 390 x 465	main unit : 300 x 360 x 530 power control unit : 330 x 160 x 125 control panel : 155 x 200 x 100
	mass ( kg )	main unit : 30 power control unit : 31	main unit : 13 power control unit : 1 control panel : 0.8



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

**JUN 21 2000**

Ms. Sheila Driscoll  
Senior Product Safety Engineer  
Quality Management Department  
Canon U.S.A., Inc.  
One Canon Plaza  
Lake Success, NY 11042-1198

Re: K992606  
Trade Name: Canon Laser Blood Flowmeter (CLBF) Model 100  
Regulatory Class: II  
Product Code: 86 HLI  
Dated: April 19, 2000  
Received: April 20, 2000

Dear Ms. Driscoll:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Pre-market Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Ms. Sheila Driscoll

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-6413. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Nancy C. Brogdon  
Acting Director  
Division of Ophthalmic and Ear,  
Nose and Throat Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Indications Statement

Page 1 of 1

510(K)Number(if known): K992606

Device Name: CLBF MODEL 100

Indications for Use:

Canon Laser Blood Flowmeter CLBF model 100 is intended for use in observing patient's fundus (retina) with the information of blood flow (velocity) in a retinal vessel. The information of blood flow is also obtained by CLBF model 100, and it is intended to help doctors' reading of retinal images.

(PLEASE DO NOT WRITE BELOW THIS LINE--CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation(ODE)

Prescription Use  (Per 21 CFR 801.109)

OR

Over-The-Counter Use

(Optional Format 1-2-96)

Denis L. Mc Carthy

(Division Sign-Off)  
Division of Ophthalmic Devices

510(k) Number K992606