

DEC 22 1999

## 510(k) Summary

---

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) number is: K9941.45

---

- 1. Submitter name, address, contact** Ortho-Clinical Diagnostics, Inc.  
100 Indigo Creek Drive  
Rochester, New York 14626-5101  
(716) 453-4041

Contact Person: Marlene A. Shulman
- 2. Preparation date** Date Special 510(k) prepared: December 7, 1999
- 3. Device name** Trade or Proprietary Name:  
VITROS Immunodiagnostic Products Free T3 Reagent Pack  
VITROS Immunodiagnostic Products Free T3 Calibrators

Common Name : Free T3 assay  
Classification Name: Total triiodothyronine test system ( 21 CFR 862.1710).
- 4. Predicate device** The VITROS Immunodiagnostic Products Free T3 Reagent Pack (new formulation) and VITROS Immunodiagnostic Products Free T3 Calibrators (new formulation) are substantially equivalent to the VITROS Immunodiagnostic Products Free T3 Reagent Pack (original formulation) and VITROS Immunodiagnostic Products Free T3 Calibrators (original formulation).

---

*Continued on next page*

## 510(k) Summary, Continued

---

- 5. Device description** The VITROS Immunodiagnostic System uses luminescence as the signal in the quantitative and semi-quantitative determination of selected analytes in human body fluids, commonly serum and plasma. Coated microwells are used as the solid phase separation system.

The system is comprised of three main elements:

1. The VITROS Immunodiagnostic Products range of immunoassay products (in this case VITROS Immunodiagnostic Products Free T3 Reagent Pack, VITROS Immunodiagnostic Products Free T3 Calibrators, which are combined by the VITROS Immunodiagnostic System to perform the VITROS Free T3 assay.
2. The VITROS Immunodiagnostic System – instrumentation, which provides automated use of the immunoassay kits. The VITROS Immunodiagnostic System was cleared for market by a separate 510(k) pre-market notification (K962919).
3. Common reagents used by the VITROS System in each assay. The VITROS Immunodiagnostic Products Signal Reagent and VITROS Immunodiagnostic Products Universal Wash Reagent were cleared as part of the VITROS Immunodiagnostic Products Total T3 Reagent Pack and VITROS Immunodiagnostic Products Total T3 Calibrators 510(k) premarket notification (K964310).

The VITROS System and common reagents are dedicated specifically for use only with the VITROS Immunodiagnostic Products range of immunoassay products.

- 6. Device intended use** VITROS Free T3 Reagent Pack  
For in vitro diagnostic use only.  
The VITROS Free T3 Reagent Pack quantitatively measures Free Triiodothyronine ( Free T3) concentration in serum and plasma (EDTA or heparin).
- VITROS Free T3 Calibrators - For *in vitro* use in the calibration of the VITROS Immunodiagnostic System for the quantitative measurement of Free T3 in serum and plasma ( EDTA or heparin).

---

*Continued on next page*

## 510(k) Summary, Continued

7. **Comparison to predicate device** The VITROS Immunodiagnostic Products Free T3 Reagent Pack (new formulation) and VITROS Immunodiagnostic Products Free T3 Calibrators (new formulation) are substantially equivalent to VITROS Immunodiagnostic Products Free T3 Reagent Pack for use with human serum which was cleared by the FDA (K970016) for in vitro diagnostic use.

Table 1 lists the characteristics of the assays performed using the VITROS Free T3 assay (new formulation) and the VITROS Free T3 assay (original formulation).

**Table 1 List of Assay Characteristics: Comparison to Predicate Device**

<b>Device Characteristic</b>	<b>New Device VITROS Free T3 assay (New formulation)</b>	<b>Predicate Device VITROS Free T3 assay (Original formulation)</b>
Sample type	Serum and plasma (EDTA or heparin).	Serum.
Antibody	One sheep polyclonal anti-T3 antibody in conjugate reagent. Immunogen: T3-edestin (a plant seed globulin). Source: Purchased from a supplier.	One sheep polyclonal anti-T3 antibody in conjugate reagent. Immunogen: T3-bovine serum albumin. Source: In-house.
Calibration range	0-35 pmol/L	0-35 pmol/L
Basic principle	Solid phase immunoassay	Solid phase immunoassay
Tracer	Enzyme labeled	Enzyme labeled
Instrumentation	VITROS Immunodiagnostic System	VITROS Immunodiagnostic System
Sample volume	25µL	25µL
Incubation time and temperature	16 minutes at 37°C with shaking	16 minutes at 37°C with shaking

*Continued on next page*

## 510(k) Summary, Continued

---

- 8. Conclusions** The information presented in the pre-market notification demonstrate that the performance of the VITROS Free T3 assay (new formulation) for use with human serum and plasma (EDTA or heparin) is substantially equivalent to the cleared predicate device.

Equivalence was demonstrated using manufactured reagents along with patient samples with measured Free T3 values spanning the assay range.

The information presented in the premarket notification provide a reasonable assurance that the VITROS Free T3 assay (new formulation) for use with human serum and plasma (EDTA or heparin) is safe and effective for the stated intended use.

---



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

DEC 22 1999

Ms. Marlene A. Shulman  
Regulatory Affairs Associate  
Ortho-Clinical Diagnostics  
100 Indigo Creek Drive  
Rochester, New York 14626

Re: K994145  
Trade Name: VITROS Immunodiagnostic Products Free T3 Reagent Pack  
VITROS Immunodiagnostic Products Free T3 Calibrators  
Regulatory Class: II  
Product Code: JIS, CPD  
Dated: December 7, 1999  
Received: December 8, 1999

Dear Ms. Shulman:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Pre-market Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895.

A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.


Page 2

Under the Clinical Laboratory Improvement Amendments of 1988 (CLIA-88), this device may require a CLIA complexity categorization. To determine if it does, you should contact the Centers for Disease Control and Prevention (CDC) at (770) 488-7655.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification"(21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597, or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Steven I. Gutman, M.D, M.B.A.  
Director  
Division of Clinical  
Laboratory Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

# Statement of Intended Use

510(k) Number (if known): K994145

Device Name: VITROS Immunodiagnostic Products Free T3 Reagent Pack  
VITROS Immunodiagnostic Products Free T3 Calibrators

Indications for Use: VITROS Free T3 Reagent Pack  
For *in vitro* diagnostic use only.  
The VITROS Free T3 Reagent Pack quantitatively measures free triiodothyronine (FT3) concentration in serum and plasma (EDTA or heparin) to aid in the differential diagnosis of thyroid disease.

VITROS Free T3 Calibrators - For *in vitro* use in the calibration of the VITROS Immunodiagnostic System for the quantitative measurement of Free T3 in serum and plasma (EDTA or heparin)

Debra Cooper  
(Division Sign-Off)  
Division of Clinical Laboratory Devices  
510(k) Number K994145

\_\_\_\_\_  
(Div  
Div  
510

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use  OR Over-The-Counter Use   
(Per 21 CFR 801.109)

Ortho-Clinical Diagnostics VITROS Immunodiagnostic Products Free T3 Reagent Pack  
VITROS Immunodiagnostic Products Free T3 Calibrators