


STAARVISC™ II
Sodium Hyaluronate

PRODUCT INFORMATION

DESCRIPTION

STAARVISC™ II is a sterile nonpyrogenic, transparent viscoelastic solution of highly purified sodium hyaluronate. **STAARVISC II** contains 12 mg/ml of high molecular weight (greater than 1 million daltons) sodium hyaluronate dissolved in physiological saline. In the limit of zero shear rate, the viscosity is 105,000 cps (105 Pa s) at 25°C, and the osmolality is approximately 320 milliosmoles (1).



Sodium hyaluronate is a naturally occurring polysaccharide composed of sodium glucuronate and N-acetylglucosamine found throughout the tissues of both man and animals. **STAARVISC II** is prepared from the dermis of rooster combs (2).

INDICATIONS

STAARVISC II is intended for use during surgery in the anterior (3) and posterior (4) segments of the human eye. Procedures include:

- Cataract extraction
- Intraocular lens (IOL) implantation
- Corneal transplantation surgery
- Glaucoma filtering surgery
- Surgical procedures to reattach the retina

STAARVISC II is designed to create and maintain anterior chamber depth and visibility, protect corneal endothelial cells and other intraocular tissues, minimize interaction between tissues during surgical manipulation, and act as a vitreous substitute during retinal reattachment surgery (5). **STAARVISC II** also preserves tissue integrity and good visibility when used to fill the anterior and posterior segments of the eye following open sky procedures.

CONTRAINDICATIONS

At the present time there are no contraindications to the use of **STAARVISC II** when used as recommended.

PRECAUTIONS

Those precautions normally considered during anterior segment and retinal attachment procedures are recommended.

Transient increases in intraocular pressure may occur following surgery because of preexisting glaucoma or due to the surgery itself (6). For these reasons, the following precautions should be considered.

- An excess quantity of **STAARVISC II** should not be used.
- **STAARVISC II** should be thoroughly removed from the anterior chamber after surgery to prevent or minimize post-operative intraocular pressure increases (spikes).

- If the postoperative intraocular pressure increases above expected values, appropriate therapy should be initiated.
- **STAARVISC II** is prepared from a biological source and the physician should be aware of the possible effects of using any biological material.
- A single use disposable cannula, such as the one provided in this package, should be used when administering **STAARVISC II**. Reuse of cannula should be avoided. The repeated use of a cannula could release particulate matter as **STAARVISC II** is injected.
- There have been isolated reports of diffuse particulates or haziness appearing after injection of products similar to **STAARVISC II** into the eye. While such reports are infrequent and seldom associated with any effects on ocular tissues, the physician should be aware of the occurrence. If observed, the particulate matter should be removed by irrigation and or aspiration.

APPLICATIONS

Cataract surgery and IOL Implantation

The required amount of **STAARVISC II** is slowly infused through a needle or cannula into the anterior chamber. The protective effect of **STAARVISC II** as a surgical aid is optimized when the injection is performed prior to cataract extraction and insertion of the IOL, and is effective for both intra- and extracapsular cataract procedures. **STAARVISC II** may be applied to the IOL prior to insertion.

Corneal transplant surgery

The corneal button is removed and the anterior chamber filled with **STAARVISC II** until it is level with the surface of the cornea. The donor graft is then placed on top of the **STAARVISC II** and sutured into place.

Glaucoma filtration surgery

STAARVISC II is injected through a corneal paracentesis to restore and maintain anterior chamber volume during the performance of the trabeculectomy.

Intraocular injection with scleral buckling procedures for retina reattachment

After release of subretinal fluid and development of buckling by tying the mattress sutures, air is injected into the vitreous cavity and then exchanged with **STAARVISC II** injected through a needle (22 to 30 gauge) passed via the pars plana epithelium. The volume of **STAARVISC II** injected (2-4ml) will vary with the volume of the subretinal fluid released and the space occupied by the buckle.

HOW SUPPLIED

STAARVISC II is a sterile viscoelastic preparation supplied in a disposable glass syringe delivering 0.8 mL of sodium hyaluronate dissolved in physiological saline. Each mL contains 12 mg of sodium hyaluronate, 9 mg of sodium chloride and q.s. sterile water for injection USP. Sodium hydroxide and/or hydrochloric acid are used to adjust pH (if necessary). **STAARVISC II** is sterile filtered and aseptically transferred to syringes. The filled syringes are sealed and the final packaging is sterilized by irradiation. Contents are sterile and nonpyrogenic in unopened and undamaged pouches. Do not use if the package is opened or damaged.

STORAGE AND HANDLING

Store at 2 to 8°C. Protect from freezing. Refrigerated **STAARVISC II** should be allowed to reach room temperature (approximately 20 to 45 minutes) prior to use.

Caution: Federal Law (USA) restricts this device to sale by, or on the order of, a physician.

REFERENCES


1. Arshinoff, Steve A., MD; "Dispersive and Cohesive Viscoelastic Materials in Phacoemulsification;" *Ophthalmic Practice*; 13:3, 1995
2. Swann DA. Studies of Hyaluronic Acid. I. The preparation and properties of rooster comb hyaluronic acid. *Biochim Biophys Acta* 1968; 156:17.
3. Pape LG, Balazs EA. The use of sodium hyaluronate (Healon®) in human anterior segment surgery. *Ophthalmol* 1980; 87:699.
4. Pruett RC, Schepens CL, Swann DA. Hyaluronic acid vitreous substitute. A six year clinical evaluation. *Arch Ophthalmol* 1979; 97:2325.
5. Miller D, Stegmann R. Use of Na-hyaluronate in anterior segment eye surgery. *AM Intra-Ocular Implant Soc J* 1980; 6:13.
6. Miller D, Stegmann R. The use of Healon® in intraocular lens implantation. *Int Ophthalmol Clinics* 1982; 22:177.

Manufactured by:
Anika Therapeutics, Inc.
Woburn, MA 01801 USA

Distributed by:
STAAR Surgical Company
Monrovia, CA 91016 USA
(626) 303-7902

P/N 490-003 rev 03/01
STAARVISC™ is a trademark of STAAR Surgical Company
© Copyright 2000 Anika Therapeutics, Inc. All Rights Reserved.

Staarvisc II sodium hyaluronate



Staarvisc II
sodium hyaluronate

P/N 490-004

Staarvisc II
sodium hyaluronate

0.8 mL
12 mg/mL

STERILE



CAUTION: For sale by or at the order of a physician.

0.8 mL
12 mg/mL




ANIKA
ANIKA THERAPEUTICS, INC.
WOBURN, MA 01801 U.S.A.

Distributed
By:


STAAR Surgical Company
Morrovia, CA 91016 U.S.A.






Starrvise II

 sodium hyaluronate





 For Single Use Only

Staar

Surgical

 P/N 490-002

CONTENTS STERILE

 IF PACKAGE NOT OPENED OR

 DAMAGED

 0.8 mL


 12 mg/mL

Manufactured By:

 Anika Therapeutics, Inc.

 Woburn, MA 01801

 U.S.A.

0.8 mL
(2 mg/mL (w)) **STERILE**
Staarvisc II
Sodium Hyaluronate
STAAR Surgical Company
Monrovia, CA 91016 U.S.A. 

(1) 27GA
Anterior Chamber Cannula

Contents sterile if pouch
is not opened or damaged.
Do not resterilize.

SINGLE USE ONLY 

STERILE | R

LOT

Manufactured For:

Anika Therapeutics, Inc.
Woburn, MA 01801 USA

P/N 300-019