

Chapter 1 Introduction

Device Description

The device is a high frequency oscillator ventilator. An earlier version of the device (3100A) has been approved for treatment of respiratory failure in infants and children. The version of the device now under review (3100B) has increased power capability and other modifications to allow treatment of adults. In a simplified description the patient circuit of the ventilator is a high-flow CPAP system. Oscillations are superimposed on the gas in the patient circuit using an electrically-driven diaphragm, similar to an audio loudspeaker cone. The oscillation frequency and magnitude can be varied. The frequency can be set between 3 and 15 cycles per second. The mean airway pressure can be set from approximately 5 to 55 cmH₂O and the bias flow (continuous sweep flow through the circuit) can be set from 0 to 60 liters per minute. The maximum pressure swing is approximately 140 cmH₂O measured at the patient circuit. Corresponding pressure swings in the trachea would be in the range of 10% of this value because of attenuation in the tracheal tube. The maximum tidal volume will be approximately 250 ml depending on the ventilator settings, tracheal tube size and the patient's pulmonary compliance. Typical settings are considerably less than these maximum values. The tidal volumes typically used are similar to the volume of the anatomic deadspace. Various mechanisms have been described to explain how these small volumes cause effective gas exchange (summarized in Krishnan and Brower, 2000).

The SensorMedics 3100 B includes alarms for overpressure and low pressures that will detect certain problems such as circuit disconnects, and some partial obstructions. The air-oxygen blender, oxygen monitor, and humidifier are connected before the gas inlet to the patient circuit; these elements are conventional and are provided by the user.

Indications for Use

The SensorMedics 3100B is indicated for use in the ventilatory support and treatment of selected patients 35 kilograms and greater with acute respiratory failure.

Contraindications

The SensorMedics 3100B Oscillatory Ventilator has no specific contraindications.

Warnings

The following Warnings must be read and understood before an attempt is made to operate the Model 3100B HFOV:

- The Model 3100B was not studied for use in children. Similar devices, SensorMedics 3100 and 3100A, are indicated for use in infant and children.
- Do not attempt to defeat the proper connection of the ground wire as it may cause damage to the device or interconnected equipment and could be injurious to the patient or to those associated with the device use. This device is factory equipped with a hospital-grade AC

power plug. Grounding reliability can only be assured when connected to a tested receptacle labeled "Hospital Grade."

- Do not operate radio transmitters within 20 feet of this instrument. This may result in erroneous pressure readings leading to false alarms and automatic shutdown.
- Do not shorten the 30" bias flow tube provided with the patient circuit as this may reduce the maximum ΔP by allowing the oscillatory pressures to be attenuated by closer proximity to the volume of the humidifier canister.
- Do not attempt to substitute a circuit configuration from any other instrument. Use of a non-3100B circuit can result in injury to the patient or to the operator, and it may cause damage to the equipment. The Patient Circuit described in this manual is specifically designed for patient use with the Model 3100B HFOV.
- There is no data to suggest that aerosols can be effectively delivered during high frequency oscillatory ventilation. Use of conventional aerosol therapy with the 3100B will probably be ineffective. Therefore, alternative vehicles for drug delivery should be considered for patients requiring this therapy.
- Only SensorMedics approved lubricants should be used. Use of any other lubricants could result in damage to the Driver Diaphragm or Bellows Water Trap Membrane causing ventilator failure or patient injury.
- The operational verification and startup procedure (Chapter 7) must be followed before ventilation of a patient commences. If at any time during the operational verification and startup procedure any abnormal function of the Model 3100B HFOV is noted, do not proceed with patient ventilation as this could cause patient injury or death; immediately contact SensorMedics Technical Support before proceeding any further.
- An audible alarm indicates the existence of a condition potentially harmful to the patient and should not go unattended. Failure to respond to alarms could result in injury (including death) to the patient and/or damage to the ventilator.
- Ensure that the cooling fan at the rear of the driver enclosure is operational.
- Under no circumstances should the ventilator be used in the presence of flammable anesthetics due to the possibility of explosion.
- Under no circumstances should a proximal airway gas temperature of 41°C be exceeded. This could result in injury to the patient's airway membranes.
- **Do Not** use the 3100B ventilator in environments where the ambient temperature is at or above 84°F (28°C). Use of the ventilator in these environments will result in extreme reduction in relative humidity in the patient's airway and possible desiccation of the patient airways.
- Failure to comply with the recommended maintenance procedures described in Chapter 8 could result in injury to the patient or operator or could result in damage to the equipment.
- Severe COPD and asthma were exclusion criteria from the randomized controlled trial of the 3100B. The benefits and/or risks associated with use of the 3100B in these patients are unknown. High frequency oscillatory ventilation is known to be less effective in diseases with increased airway resistance and its use may potentially result in air trapping and hyperinflation. This should be taken into consideration if used in these patients.

Precautions

The following Cautions must be read and understood before an attempt is made to operate the Model 3100B HFOV:

- Follow closely the recommendations contained in Chapter 9, Clinical Guidelines, regarding the use of chest radiographs to monitor patient condition. During HFOV, as with all ventilators, the relationship between improvement in lung compliance, inadvertent increases in lung volume, increased pleural pressure, and decreased venous return is a matter of concern, since it may result in decreased cardiac output.
- Patient size is an important guideline as to lung volume and anatomical dead space, as well as the metabolic demand placed on ventilation. While the maximum displacement volume of the 3100B is approximately 365 ml, the actual volume delivered to the patient is dependent on power setting, frequency, endotracheal tube size, and patient respiratory system compliance. It is recommended that the operator review Chapter 9 of this manual, "Clinical Guidelines."
- The patient's $t\text{PCO}_2$ and $t\text{PO}_2$ or SpO_2 should be monitored continuously to insure that blood gases are at the proper level. It is important that an unrestricted and unobstructed patient airway be maintained during HFOV. To insure a patent airway, always maintain proper suctioning procedures as described in the Suctioning Guidelines Section of Chapter 9, Clinical Guidelines. Since only proximal airway pressure is measured, no alarm will occur in the event of an obstruction or restriction.
- Ensure that the stopcock is closed prior to performing a Patient Circuit Calibration. If the Water Trap Stopcock is left open, Patient Circuit Calibration (39-43 cmH₂O) may not be achievable, and the deliverable Pa will be reduced.
- Deviation from the assembly methods described in Chapter 6, Assembly and Installation, could damage the Model 3100B, render it mechanically unstable, or cause it to malfunction. If any questions arise regarding the assembly procedure, please contact SensorMedics Technical Support immediately before proceeding.
- Care should be taken not to crimp or perforate any of the control or sensing lines (running to or from the Patient Circuit) during assembly, operating or cleaning of the ventilator as this will cause malfunction of the Safety Alarms, Warning Alarms, Caution Alarms, and/or Pressure Limit Controls.
- Before attaching the patient circuit to the ventilator, the driver diaphragm of the 3100B should be inspected for cuts and tears. If any damage is noted, do not proceed with patient ventilation as this could cause failure of the ventilator. Immediately contact SensorMedics Technical Support for assistance.
- The driver diaphragm of the 3100B has been coated with a special lubricant during assembly. Please do not clean the driver diaphragm with cleaning solvents as it may degrade the materials causing premature wear of the driver diaphragm.
- When connecting the Patient Circuit, make certain that it is properly supported and oriented by the support arm as described in Chapter 6, Assembly and Installation. Failure to do so could result in inadvertent Patient Circuit disconnection due to oscillatory forces or could result in collection of humidifier condensate in the patient airway.
- If the temperature probe is wiped with alcohol, allow the alcohol to evaporate completely before inserting it into the circuit. A high residual of alcohol can weaken the acrylic adapter and cause fracturing.

- Proper operation of the ventilator must be verified prior to each use. Refer to Chapter 7, Operational Verification and Startup Procedures. The alarm functions tested in this procedure verify the capability of the device to detect and indicate conditions which could have a harmful effect on the patient.
- Touch the outer metal cabinet of the instrument before touching any other component to avoid possible instrument component damage from Electrostatic Discharge.
- When the ventilator is connected to a patient, it is imperative that someone be in attendance at all times in order to react to any alarms and to detect other indications of a problem.
- The Inlet Filter Cartridges for the blended gas and the air inputs to the ventilator must be changed at least every 500 hours of operation as described in Chapter 8, Maintenance and Troubleshooting. Failure to replace a Filter Cartridge or substitution of an unauthorized cartridge could result in injury to the patient and/or damage to the equipment. Use only SensorMedics Inlet Filter Cartridges.
- The filter cartridge body must be screwed back on securely. Cross-threaded or loose installation will result in leaks and possible dislodging of the cartridge body. If the cartridge body is dislodged, it will cause the ventilator to cease functioning.
- The cover enclosing the Control Package, Column, or any other portion of the ventilator must not be removed by the user. To avoid electrical shock hazard, please refer all service requiring cover removal to a qualified biomedical equipment service technician.
- Recheck and readjust alarm levels after any parameter change has been made.
- Troubleshooting with the 3100B should be done "OFF PATIENT" to avoid any potentially dangerous situations such as abrupt changes in the Pa.
- Do not use extraneous ventilator circuit attachments (such as a suction port) without a secondary external alarm capable of detecting ventilator disconnection. Due to their inline pressure characteristics, such attachments could possibly keep the Pa alarm from detecting an accidental ventilator circuit disconnection.
- Fractional concentration of inspired oxygen should be verified with an oxygen monitor. Administration of excessive oxygen to a patient may be harmful. It is imperative that the prescribed gas mixture is delivered by the blending system.
- The Water Trap must be drained at intervals as described in Chapter 8, Maintenance and Troubleshooting. If the ventilator is operating, leave a small amount of water at the bottom of the Water Trap container to act as a flow and pressure seal between the ventilator and the output of the drain.
- To help prevent patient injury due to humidifier malfunction, the use of a humidifier with the following characteristics is strongly recommended:
 - Thermally protected heater.
 - Alarms on over-filled water reservoir.
 - Alarms on under-filled water reservoir.
 - Alarms when electrically open or shorted temperature probe detected.
 - Alarms at probe temperatures > 41°C.
 - Alarms when dislodged temperature probe detected.
- Do not place on the Control Package of the ventilator any fluid-containing accessories, accessories that weigh more than ten pounds, or accessories that extend more than six inches

above the ventilator electronics package or beyond its sides. This could cause damage to the ventilator, or could cause the ventilator to tip over, resulting in patient or user injuries and/or damage to the equipment.

- Do not overturn the Patient Circuit Calibration screw, as this may cause damage to the device. When it is nearing its adjustment limit, it will reach a mechanical stop.
- Do not allow liquids to penetrate the air vents of the ventilator as this may result in machine failure or malfunction.
- Do not use a liquid sterilization agent on the outside of the ventilator as this may cause damage.

Adverse Effects

Observed Adverse Events

A prospective clinical study of the Model 3100B in patients with acute respiratory distress syndrome (ARDS) was conducted at ten sites. The 148 patients enrolled in the study were randomly assigned to one of two groups: a treatment (i.e., "high-frequency") group, in which patients were treated with the Model 3100B; or a control (i.e., "conventional") group, in which patients were treated with a conventional ventilator. Treatment outcomes and adverse events were determined at one and six months. Table 1 summarizes the adverse events observed during the study.

Table 1: Observed adverse events







Adverse event	# (%) of patient who had this event in the treatment group, N=75
Mucus-plugged ET tube	2 (3%)
Inadequate Oxygenation	4 (5%)
Respiratory acidosis	4 (5%)
New or worsening air leak syndrome	7 (9%)
Intractable hypotension	0 (0%)

Potential adverse events

The adverse events associated with the use of high-frequency ventilation include: atelectasis, inadequate oxygenation, intractable hypotension, mucus-plugging of the tracheal tube, necrotizing tracheobronchitis, new or worsening air leak syndrome, over-humidification, under-humidification, over-ventilation, under-ventilation, pneumothorax, pneumopericardium, pneumomediastinum, pneumoperitoneum, pulmonary interstitial emphysema and respiratory acidosis.

Symbols

The following symbols are used on this device:

Symbol	Compliance	Meaning
	Symbol #03-02 IEC 60878	Indicates ATTENTION, consult ACCOMPANYING DOCUMENTS
	Symbol # 5333 IEC 60417 Symbol #02-03 IEC 60878	This symbol indicates TYPE B equipment, which indicates equipment that provides a particular degree of protection against electric shock, particularly with regards to allowable leakage current and reliability of the protective earth connection.
	Symbol #5032 IEC 60417 Symbol #01-14 IEC 30878	This symbol indicates the equipment is suitable for alternating current.
	Symbol #5007 IEC 60417 Symbol #01-01 IEC 60878	Indicates circuit breaker ON (Power)
	Symbol #5008 IEC 60417 Symbol #01-02 IEC 60878	Indicates circuit breaker OFF (Power)
	SensorMedics Symbol	Position Lock. Clockwise rotation locks instrument top. Counter-clockwise rotation unlocks instrument top, allowing it to be swiveled for best view of front controls and displays.

Exterior Labels

This section identifies the labels attached to the exterior of the 3100B. All labels are shown at approximately their actual size. Your system may not have all of the labels listed.

Patient Circuit Calibration





Figure 1.1. Patient Circuit Calibration Procedure Label.

The Patient Circuit Calibration Procedure Label describes the steps necessary to calibrate the patient circuit to the 3100B. This procedure is also explained in the Patient Circuit Calibration Section of Chapter 8, Maintenance and Troubleshooting.

Ventilator Performance Checks



767165-101F

Figure 1.2. Ventilator Performance Checks Label.

The Ventilator Performance Checks Label assists in setting Power, Mean Pressure Adjust, and Bias flow controls to achieve specific ranges of ΔP and P_a . These procedures are explained in the Performance Verification Section of Chapter 7, Operational Verification and Start-up Procedures.

Blender/Cooling Gas Filter Replacement Record

BLENDER COOLING GAS FILTER REPLACEMENT RECORD			
<i>Filters should be replaced every 500 hours of operation.</i>			
REPLACED	REPLACED	REPLACED	REPLACED
____ Hrs	____ Hrs	____ Hrs	____ Hrs
____ Hrs	____ Hrs	____ Hrs	____ Hrs
____ Hrs	____ Hrs	____ Hrs	____ Hrs
____ Hrs	____ Hrs	____ Hrs	____ Hrs
____ Hrs	____ Hrs	____ Hrs	____ Hrs
____ Hrs	____ Hrs	____ Hrs	____ Hrs
____ Hrs	____ Hrs	____ Hrs	____ Hrs
____ Hrs	____ Hrs	____ Hrs	____ Hrs
____ Hrs	____ Hrs	____ Hrs	____ Hrs
____ Hrs	____ Hrs	____ Hrs	____ Hrs
____ Hrs	____ Hrs	____ Hrs	____ Hrs
____ Hrs	____ Hrs	____ Hrs	____ Hrs
____ Hrs	____ Hrs	____ Hrs	____ Hrs
____ Hrs	____ Hrs	____ Hrs	____ Hrs

P/N 767168 A

Figure 1.3. Blender/Cooling Gas Filter Replacement Record Label.

The Blender/Cooling Gas Filter Replacement Record Label provides a place to document the 500 hour gas filter changes. For more information, see the Operator Maintenance Procedures Section of Chapter 8, Maintenance and Troubleshooting.

Driver Replacement Record

Driver Replacement Record			
<i>Drivers should be replaced every 4000 hours of operation.</i>			
Replaced	Driver S/N	Replaced By	Date
____ Hrs			
____ Hrs			
____ Hrs			
____ Hrs			
____ Hrs			
____ Hrs			
____ Hrs			
____ Hrs			
____ Hrs			
____ Hrs			
____ Hrs			
____ Hrs			
____ Hrs			
____ Hrs			
____ Hrs			

P/N 768338A

Figure 1.4. Driver Replacement Record Label.

The Driver Replacement Record Label provides a place to document the 4,000 hours replacement of the Oscillator Subassembly. For more information, see the Scheduled Periodic Maintenance Section of Chapter 8, Maintenance and Troubleshooting.

Radio Frequency Interference (RFI) Warning

– WARNING –

DO NOT OPERATE RADIO-TRANSMITTERS WITHIN 20 FEET OF THIS INSTRUMENT. THIS MAY RESULT IN ERRONEOUS PRESSURE READINGS LEADING TO FALSE ALARMS AND AUTOMATIC SHUT-DOWN.

– SEE OPERATOR'S MANUAL –

P/N 768559C

Figure 1.5. Radio Frequency Interference Warning Label.

The Radio Frequency Interference (RFI) Warning Label refers to the possible problems caused by interference from hand-held radio transmitters. The RFI warning is also discussed in the Troubleshooting Section of Chapter 8, Maintenance and Troubleshooting.

Name Rating Label



3100B OSCILLATORY VENTILATOR	
VOLTAGE	115V  HERTZ 60
AMPERE	7.5A
S/N	31003
CAT P/N	770155
	SENSORMEDICS MADE IN U.S.A.
	SensorMedics Corporation 22705 Savi Ranch Parkway Yorba Linda, California 92687

Figure 1.6. Name Rating Label.

The Name Rating Label lists specific information on each individual instrument: the Model Name and Number, the Voltage and Current Rating, the Serial Number, and the instrument's Catalog Part Number. (The example shown is for the 115V, 7.5A, 60Hz model; your instrument may have a different rating.)

Battery Attachment




Figure 1.7. Battery Attachment Label.

The Battery Attachment Label indicates the correct position for the installed Power Failure Alarm Battery. For directions on changing the battery, see the Changing the Power Failure Alarm Battery Section of Chapter 8, Maintenance and Troubleshooting.

Battery Specification



Figure 1.8. Battery Specification Label.

The Battery Specification Label indicates the type of Power Failure Alarm Battery (9V alkaline) that must be used. For directions on changing the battery, see the Changing the Power Failure Alarm Battery Section of Chapter 8, Maintenance and Troubleshooting.

Chapter 2 Clinical Study

A prospective clinical study of the Model 3100B in patients with acute respiratory distress syndrome (ARDS) was conducted at ten sites. The 148 patients enrolled in the study were randomly assigned to one of two groups: a treatment (i.e., "high-frequency") group, in which patients were treated with the Model 3100B; or a control (i.e., "conventional") group, in which patients were treated with a conventional ventilator.

Inclusion and Exclusion Criteria

Patients were eligible for inclusion in the study if they met the following criteria:

- at least 16 years of age;
- at least 35 kilograms;
- $PaO_2/FiO_2 < 200$;
- bilateral pulmonary infiltrates not resulting from left atrial hypertension; and
- positive end-expiratory pressure (PEEP) of at least 10 cm H₂O.

Patients were excluded if any of the following were true:

- informed consent could not be obtained;
- patient had been treated with FiO_2 greater than 80% for at least 48 hours;
- patient had severe persistent air leak;
- patient had non-pulmonary terminal prognosis;
- patient had severe chronic obstructive pulmonary disease
- patient had asthma; or
- patient had been recently enrolled in another ARDS or septic shock investigation.

Methods

The general treatment goal for the high-frequency group and the conventional group was the same: to maintain an O₂ saturation of at least 88%; and to maintain a pH of greater than 7.15, while minimizing peak pressures and treating metabolic acidosis. The mean airway pressure was maintained until FiO_2 had been reduced to less than 60%, after which mean airway pressure and FiO_2 were given equal priority for reduction.

In the high-frequency group, the Model 3100B was initially set to provide pressure oscillations at a frequency of 5 Hz, with the mean airway pressure set 5 cm H₂O higher than the ventilator setting used before the patient was enrolled in the study, and with the oscillation amplitude (ΔP) set for adequate chest wall vibration. If ventilation was inadequate, ΔP was increased. If ventilation was still inadequate with maximum ΔP , the frequency of pressure oscillations was reduced in 1 Hz steps. When mean airway pressure had been decreased to less than 30 cm H₂O, or when there was no progress in weaning with the Model 3100B, ventilator weaning continued using a conventional ventilator. A patient assigned to the high-frequency group was treated with the Model 3100B protocol until: consent was withdrawn; the patient had died or been weaned from mechanical ventilation; the patient had been ventilated for 30 days; or the patient met defined treatment failure criteria and would, in the opinion of their physician, benefit from conventional ventilation.

In the conventional group, patients were treated with conventional pressure-control, volume-limited ventilation with an inspiratory to expiratory (I:E) ratio of approximately 1:2. Tidal volumes nominally between 6 and 10 mL/kg were used. (The average tidal volume delivered was 10.2 mL/kg of ideal body weight.) If oxygenation was inadequate, PEEP was increased in increments of up to 5 cm H₂O to improve oxygenation. If oxygenation was still inadequate with PEEP greater than or equal to 18 cm H₂O, the I:E ratio was increased incrementally. A patient assigned to the conventional ventilation was treated within the study with the conventional ventilator until: consent was withdrawn; the patient had died or been weaned from mechanical ventilation; or the patient had been ventilated for 30 days.

Patient outcomes were determined after one month and six months. The possible outcomes were:

- death;
- survival with respiratory support; or
- survival without respiratory support.

In this trial, "respiratory support" was defined to include mechanical ventilation, CPAP or supplemental oxygen. Only survival without respiratory support was considered a successful outcome. Outcome data were analyzed using an "intention to treat" analysis.

Hypothesis

The primary hypothesis was that the proportions of patients in the high-frequency and conventional groups with unsuccessful one-month outcomes would be equivalent.

Statistically stated, the hypothesis was that the proportion of patients with an unsuccessful one-month outcome—i.e., the proportion of patients who died or were still receiving respiratory support after one month—would be no more than 10% greater in the high-frequency group than in the conventional group, with 95% confidence.

Study Population

In this trial, 75 patients were assigned to the high-frequency group, and 73 patients were assigned to the conventional group. The patient demographics, the pre-enrollment diagnoses, the pre-enrollment ventilator settings and the pre-enrollment clinical indicators were similar for the two groups (Table 3).

Table 3: Patient demographics, pre-enrollment ventilator settings and pre-enrollment clinical indicators.

Category	Parameter	Mean \pm std. dev.		Units
		Treatment group	Control group	
Patient demographics	Age	48 \pm 17	51 \pm 18*	years
	Weight	78 \pm 25	81 \pm 26*	kg
	Gender (% male)	52	64	%
Pre-enrollment ventilator settings	Peak inspiratory pressure	39 \pm 7	38 \pm 8	cm H ₂ O
	Positive end-expiratory pressure	13 \pm 3	14 \pm 3	cm H ₂ O
	Mean airway pressure	22 \pm 5*	24 \pm 7*	cm H ₂ O
	Tidal volume per kilogram of actual body weight	8.2 \pm 3*	7.8 \pm 3*	mL/kg
	FiO ₂	71 \pm 19	72 \pm 19	%
Pre-enrollment clinical indicators	PaO ₂	76 \pm 20	73 \pm 18	mm Hg
	PaCO ₂	44 \pm 12	45 \pm 12	mm Hg
	pH	7.37 \pm 0.09*	7.34 \pm 0.11	
	PaO ₂ /FiO ₂	114 \pm 37	111 \pm 42	
	Oxygenation index ¹	24 \pm 15*	27 \pm 19*	
	Mean blood pressure	80 \pm 14*	76 \pm 12*	mm Hg
	Cardiac output	7 \pm 2*	7 \pm 3*	L/min
APACHE II score ²	22 \pm 6*	22 \pm 9*		

* Denotes that values for this parameter were not available for all patients.

¹ Oxygenation index = $100 \times \text{mean airway pressure} / (\text{PaO}_2/\text{FiO}_2)$

² APACHE II is a disease severity score.

Results

The one-month outcomes in this trial are summarized in Table 4, below. The patients in the high-frequency group had unsuccessful one-month outcomes with greater frequency than did the patients in the conventional group. Based on the 95% confidence interval computed from the one-month outcomes, the treatment in the high-frequency group could fail as much as 20% more often. Therefore, the prospectively defined hypothesis was not met. However, the mortality rate in the high-frequency group was lower than the mortality rate in the conventional group. The observed six-month outcomes are summarized in Table 5. Both the unsuccessful treatment rate and the mortality rate were lower in the high-frequency group than in the conventional group. The lower one-month and six-month mortality rates, and the lower six-month unsuccessful treatment rate, in the high-frequency group provide reasonable assurance that the Model 3100B is safe and effective.

Table 4: One-month outcomes

One-month outcome	Treatment group, N=75	Control group, N=73	Difference	
			Absolute	95% CI
Unsuccessful	78%	73%	+5%	-10% to +20%
Death	37%	52%	-15%	

Table 5: Six-month outcomes

Six-month outcome	Treatment group, N=75	Control group, N=73
Unsuccessful	47%	62%
Death	47%	59%

Observed Treatment Failures

Likely types of treatment failure were prospectively identified. The frequency with which each type of failure occurred was similar for the two groups, as shown in Table 6, below.

Table 6: Observed treatment failures

Treatment failure	Treatment group, N=75	Control group, N=73
Mucus-plugged ET tube	3%	1%
Inadequate oxygenation	5%	8%
Respiratory acidosis	5%	8%
New or worsening air leak syndrome	9%	12%
Intractable hypotension	0%	2%

Causes of Death

The causes of death were identified for those patients who died while being treated with a ventilator. Deaths due to withdrawal of mechanical ventilation were also identified. For many patients, more than one cause of death was identified. The causes of death included cardiac arrhythmia, multiple organ failure, sepsis and profound hypoxemia. The causes of death in each group occurred with similar frequency (Table 7).

**Table 7: Causes of death observed for patients
who died while being treated with a ventilator**

Causes of death	# (%) of patients who died while being treated with a ventilator	
	Treatment group, N=75	Control group, N=73
Total	8 (11%)	15 (21%)
Cardiac arrhythmia	3 (4%)	6 (8%)
Multiple organ failure	3 (4%)	5 (7%)
Sepsis	6 (8%)	8 (11%)
Profound hypoxemia	3 (4%)	6 (8%)
Other	1 (1%)	2 (3%)