

Daig Corporation

a St. Jude Medical company

LIVEWIRE TC™

Steerable Electrophysiology Catheter

INSTRUCTIONS FOR USE

Table of Contents

	<u>Page</u>
1 DEVICE DESCRIPTION	3
2 INDICATIONS AND USAGE	3
3 CONTRAINDICATIONS.....	3
4 WARNINGS AND PRECAUTIONS	4
5 ADVERSE EVENTS	6
5.1 Observed Adverse Events	6
Table 1. Observed Adverse Events	7
5.2 Potential Adverse Events	7
6 CLINICAL STUDIES.....	8
Table 2. Patient Age and Acute Success by Indication.....	9
Table 3. Recurrence-free survival (Kaplan-Meier Estimates).....	9
7 PATIENT SELECTION AND TREATMENT.....	9
7.1 Individualization of Treatment	9
7.2 Specific Patient Populations	10
8 PATIENT COUNSELING INFORMATION.....	10
9 SUGGESTED DIRECTIONS FOR USE	10
9.1 Physician Training	10
9.2 Directions for Use	10
9.3 Compatible RF Generators and Accessories.....	11
Table 4. Specifications for a Compatible RF Generator:	11
10 HOW SUPPLIED	12
10.1 Packaging.....	12
10.2 Storage and Shelf-Life.....	12
10.3 Warranty and Replacement Policy.....	12

1 DEVICE DESCRIPTION

The Daig LivewireTC™ Steerable Electrophysiology Catheter (Livewire Catheter) is a flexible electrode catheter which records electrophysiological activity from the heart and transmits radio frequency (RF) current to the catheter tip electrode for ablation purposes. The catheter is constructed of a radiopaque polyurethane shaft and platinum electrodes. All electrodes may be used for recording and stimulation, but only the tip electrodes may be used to deliver RF energy. The remote control handle located at the proximal end of the catheter deflects distal the tip of the catheter.

The Livewire Catheter must be used with a dispersive pad (reference electrode), and a compatible RF generator with the following specifications:

- A thermocouple (for temperature measurement)
- Maximum temperature limit - 100°C
- Modes - Temperature control and monitoring and power control.
- Maximum output - 50 watts
- RF output frequency - 450–550 KHz
- Impedance range - 40–250Ω

2 INDICATIONS AND USAGE

The Livewire Catheter is indicated for cardiac electrophysiological mapping and for use with a compatible RF generator for:

- interruption of accessory atrioventricular (AV) conduction pathways associated with tachycardia;
- the treatment of AV nodal re-entrant tachycardia (AVNRT); or
- creation of complete AV nodal block in patients with a difficult to control ventricular response to an atrial arrhythmia.

3 CONTRAINDICATIONS

Do not use this device:

- in patients with active systemic infection;
- via the transseptal approach in patients with left atrial thrombus or myxoma, or interatrial baffle or patch; or
- in patients with aortic valve replacement via the retrograde transaortic approach.

4 WARNINGS AND PRECAUTIONS

Warning

- **Do not ablate from within a coronary artery** as the resulting myocardial injury can be fatal. Adequate fluoroscopic visualization is necessary during the transaortic approach to avoid placement of the ablation catheter in the coronary vasculature.
- **Stroke or myocardial infarction** may occur in patients undergoing **left-sided ablation procedures**. Patients should be closely monitored during the post-ablation period for clinical manifestations of embolic events.
- **Precautions in patients with implantable pacemakers and implantable cardioverter/ defibrillators (ICDs):**
 - **Deactivate ICDs** as they could discharge and injure the patient or be damaged by the ablation procedure
 - Have temporary external sources of pacing and defibrillation available
 - **Do not apply RF energy directly to a lead** or to tissue immediately in contact with a lead because it could potentially damage the lead or lead function
 - Perform a complete analysis of the implanted device function after ablation
- **Complete AV block** can occur when ablating septal accessory pathways or in the treatment of AVNRT. Closely monitor AV conduction during RF energy delivery and immediately terminate energy delivery if partial or complete AV block is observed.
- **Minimize X-ray exposure** – Significant x-ray exposure can result in acute radiation injury as well as dose-related risk for somatic and genetic effects. Take all appropriate measures to minimize x-ray exposure to both patients and clinical staff.
- **X-ray exposure to children** – The long-term risk of protracted fluoroscopy has not been established. Therefore, careful consideration must be given for the use of the device in prepubescent children.
- **Pregnancy** – Careful consideration should be given to the use of this device in pregnant women because of the dose-related risk for somatic and genetic effects of x-ray exposure.
- **Training** – Cardiac ablation procedures should be performed only by appropriately trained personnel in a fully equipped electrophysiology laboratory (see section 10).
- **Instructions for Use** – Do not attempt to operate the Livewire Catheter before completely reading and understanding the applicable directions for use.
- **Long-term risks of RF ablation** – The long-term risks of lesions created by RF ablation have not been established. In particular, any long-term effects of lesions in proximity to the specialized conduction system or coronary vasculature are unknown.

4.1 Precautions Specific to the Livewire Catheter

- Do not use the Livewire Catheter for long term pacing.

- Use only Daig Electrophysiology connector cables with the Livewire Catheter since patient or operator injury may occur if incompatible cables are used.

4.2 Handling and Sterilization Precautions

- The Livewire is for SINGLE USE ONLY. Do not resterilize or reuse.
- Do not use the Livewire Catheter after the expiration date because the device performance may no longer be acceptable and/or the device may no longer be sterile.
- Inspect the packaging and catheter prior to use. If the package or the catheter appears damaged, do not use and contact your local Daig representative.
- Do not allow moisture to contact the connector as equipment malfunction or operator injury may occur.

4.3 Precautions During Catheter Use

- The patient should not contact grounded metal surfaces. Use only isolated amplifiers, pacing equipment, and ECG equipment or patient injury or death may occur. Leakage current from any connected device to the patient must not exceed 10 microAmps (μA) under any circumstances.
- Do not use excessive force to advance or withdraw the catheter. Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade.
- Do not insert or withdraw the catheter without straightening the catheter tip (pulling the thumb knob back)
- Do not use the catheter if the small vent area at the connector end of the handpiece is clogged since air may be forced into the catheter lumen and into the bloodstream.
- Use both fluoroscopy and electrograms to monitor the advancement the catheter to the area of the endocardium under investigation to avoid vascular or cardiac damage.

4.4 Environmental and EMI

- Catheter materials are not compatible with magnetic resonance imaging (MRI).
- During ablation procedures, this catheter is used in conjunction with an RF generator. Electromagnetic interference (EMI) produced by the radiofrequency (RF) generator during the delivery of RF power may adversely affect the performance of other equipment.

4.5 Precautions During Ablation

- Do not increase power before checking for lead connection and appropriate dispersive electrode application. Effective contact between the patient and the dispersive electrode must be verified whenever the patient is repositioned.
- Do not deliver RF energy with catheter outside the target site. The RF generator can deliver significant electrical energy and may cause patient or operator injury.

- Avoid use of electrodes and probes of monitoring and stimulating devices which could provide paths for high frequency current. Reduce the burn hazard by placing the electrodes and probes as far away as possible from the ablation site and the dispersive electrode.
- In the event of a generator cutoff (impedance or temperature), the catheter must be withdrawn and the tip electrode cleaned of coagulum before RF current is re-applied. Use only sterile saline and gauze pad to clean the tip.
- Do not scrub or twist the tip electrode as damage may cause catheter failure or patient injury.
- Discontinue ablation immediately and replace catheter if tip temperature fails to rise during ablation (temperature sensing model).
- The temperature sensing model of the catheter measures electrode tip temperature, not tissue temperature. If the generator does not display temperature (temperature sensing model), verify that the appropriate cable is plugged into the generator. If temperature still is not displayed, there may be a malfunction in the temperature sensing system which must be corrected prior to applying RF power.

5 ADVERSE EVENTS

5.1 Observed Adverse Events

The Livewire Catheter was studied in 329 patients undergoing electrophysiologic (EP) mapping and ablation in two clinical studies: 236 patients using an investigational RF generator and 93 patients using a commercially available generator. Of the 329, eleven patients were not treated with the Livewire Catheter due to either investigational generator issues (9) or non-availability of catheter accessories (2). There was 1 major (atrial perforation) and no minor adverse events in these 11 patients. These 11 patients were excluded from further adverse event analysis. The 318 patients included 10 patients treated for typical atrial flutter.

The number of patients with adverse events (major or minor) was 9 of 226 (4.0 %) in the study using the investigational generator and 2 of 92 (2.2%) in the study using the commercially available generator.

Table 1. Observed Adverse Events
 All Patients Treated (N=318), includes 10 patients with atrial flutter

Adverse Event Category Description	Total No. of Pts. Experiencing Adverse Events	% of Pts Experiencing Adverse Events [95% C.I.]
Major; Heart block (requiring a pacemaker)	4	1.3% [0.3%, 3.2%]
Minor (included puncture site bleeding, pericardial effusion, left bundle branch block, transient complete heart block, and first degree atrioventricular block)*	7	2.2% [0.9%, 4.5%]

* Patients with both minor and major adverse events (N=2) were counted only as major adverse events.
 C.I. – Confidence intervals by the exact (binomial) method

5.2 Potential Adverse Events

Adverse events (in alphabetical order) which may be associated with catheterization and ablation include:

Catheterization/catheter procedure related:

- air embolism
- arrhythmias
- AV fistula
- cardiac perforation
- hemothorax
- nerve palsy or weakness
- pneumothorax
- pseudoaneurysm
- tamponade
- thrombi
- thromboembolism
- thrombosis
- valvular damage
- vascular bleeding/local hematomas
- vasovagal reactions
- visual blurring

RF ablation related:

- cardiac perforation/tamponade
- cardiac thromboembolism
- cerebrovascular accident (CVA)
- chest pain/discomfort
- complete heart block
- coronary artery dissection
- coronary artery spasm
- coronary artery thrombosis
- pericarditis
- transient ischemic attack (TIA)
- valvular damage
- ventricular tachyarrhythmia

6 CLINICAL STUDIES

A prospective, multicenter (10 centers, 7 in the United States and 3 in Canada) study was conducted of radio frequency (RF) ablation of atrioventricular (AV) accessory pathways (AP) associated with tachycardia, AV nodal re-entrant tachycardia (AVNRT), or creation of complete AV nodal block in patients with difficult to control ventricular response to an atrial arrhythmia.

Methods: Acute success was defined as ablation of the target site and elimination of the arrhythmia with no recurrence before patient discharge. Recurrence-free survival to 12 months was assessed.

Patients Studied: The Livewire Catheter was studied with a commercially available generator (N = 93), and with an investigational RF generator, (N = 236) for patients undergoing EP mapping and RF ablation for supraventricular tachycardia (SVT). The closed loop temperature control mode was used for 92 of the patients. The procedure was aborted in 11 patients and 10 patients were found to be in typical atrial flutter. Acute effectiveness was thus assessed in 308 patients.

Results: The patients ranged in age from 11 to 84 years (mean 44 years) and included 183/329 females (56%). The acute success rate for all arrhythmias was 201 of 216 (93%) in the study using the investigational generator and 86 of 92 (93%) in the study using the commercially available generator. The difference between the results for the two generators is 0 % with a 95% confidence interval of [-9%, 8%], so the results for the two generators have been combined. Table 2 shows patient age, enrollment and acute success by indication.

Table 2. Patient Age and Acute Success by Indication
All Patients Treated for SVT (N=308)

Arrhythmia	Age, mean (range)	Success (%) [95% CI]
AVNRT	46 (11, 78)	156/163 (96%) [91%, 98%]
AP	36 (12, 66)	99/112 (88%) [82%, 94%]
AV nodal	57 (22, 84)	31/33 (94%) [80%, 99%]
All arrhythmias	44 (11, 84)	286/308 (93%) [89%, 95%]

C.I. – Confidence intervals by the exact (binomial) method

Of the 201 patients successfully treated with the investigational generator, follow-up data was available for 195. Of these 195 patients, recurrence was reported in 41 (21%), and 16 (8%) were lost to follow-up by 12 months. Table 3 shows the recurrence free survival to 12 months.

Table 3. Recurrence-free survival (Kaplan-Meier Estimates)
Patients Successfully Treated with follow-up (N=195*)

Follow-up Time	Cumulative Survival [95% CI]	Number Remaining
One Month	92% [88%, 96%]	179
Three Months	87% [82%, 91%]	164
Six Months	82% [77%, 88%]	150
Twelve Months	77% [70%, 83%]	98

C.I. – Confidence intervals from the Kaplan Meier estimates

7 PATIENT SELECTION AND TREATMENT

7.1 Individualization of Treatment

Antiplatelet or Anticoagulation Use

To avoid thromboemboli, intravenous heparin is used when entering the left heart during ablation, and many physicians prescribe aspirin, less often warfarin, for about 3 months afterward. No consensus yet exists about the need for short-term anticoagulation after ablation.

Left Heart Insertion

During the clinical study, systemic anticoagulation before intracardiac RF catheter ablation in the left heart was typically an initial intravenous heparin bolus

of 3000 - 5000 Units. Anticoagulation was maintained with an intravenous heparin drip or additional periodic intravenous boluses of heparin as necessary.

Right Heart Insertion

During the clinical study, systemic anticoagulation was variable for patients undergoing intracardiac RF catheter ablation in the right heart. If used, systemic anticoagulation in the clinical study before ablation was typically an initial intravenous heparin bolus of 2000 - 5000 Units. For some patients, anticoagulation was maintained by intravenous heparin drip or additional periodic intravenous boluses.

Choosing Temperature or Power Control Mode

Please refer to the compatible RF generator's Directions for Use for information in choosing between temperature or power control modes.

7.2 Specific Patient Populations

The safety and effectiveness of cardiac ablation has not been established in:

- Asymptomatic patients;
- patients who are pregnant; or
- nursing mothers.

8 PATIENT COUNSELING INFORMATION

Patients may require anticoagulation and/or antiplatelet therapy for an indefinite period based on the patient's condition.

9 SUGGESTED DIRECTIONS FOR USE

9.1 Physician Training

Physicians must be familiar with the techniques and appropriately trained for cardiac mapping and ablation procedures. All mapping and ablation procedures must be performed in a fully-equipped electrophysiology laboratory.

9.2 Directions for Use

1. Use a Daig Fast-Cath™ Introducer to insert the Livewire Catheter.
2. Always use fluoroscopy when positioning the electrode catheter.
3. Use only a Daig Electrophysiology ablation extension cable to connect the Livewire Catheter to the appropriate electrode interface in the electrophysiology lab.
4. To record intracardiac electrograms, connect extension cable to catheter. Observe polarity of proximal end connector pins of extension cable when connecting to an ECG amplifier of your choice.

5. To use this device for temporary pacing, connect extension cable to catheter. Observe polarity of proximal end connector pins of extension cable when connecting to an external pulse generator of your choice.
6. To use this device for RF ablation without temperature monitoring, connect extension cable to catheter. Observe polarity of proximal end connector pins of extension cable when connecting to an ECG amplifier of your choice. Determine the exact location of the ablation site using the physicians clinical experience with electrophysiological and fluoroscopic guidance and mapping techniques. Connect the proximal end connector of the extension cable to the compatible RF generator previously described in these instructions.
7. When using the Livewire Catheters with the temperature monitoring feature, a six conductor ablation extension cable (blue) must be used.
8. Consult the compatible RF generator instructions for use for the proper Connection of the patient grounding plate.
9. To manipulate the tip portion of this catheter, rotate the control collar located in the handle of the catheter.
10. Always use fluoroscopy when manipulating tip of catheter.
11. Always straighten catheter tip before removing catheter from patient.

9.3 Compatible RF Generators and Accessories

The Livewire Catheter should be used only with a RF generator which has been shown to be safe and effective for cardiac ablation.

Table 4. Specifications for a Compatible RF Generator:

PARAMETER	SPECIFICATION
Thermometry	Thermocouple
Temperature Limit, maximum	100°C
Modes: (must operate in all 3 modes)	1. Temperature Control 2. Temperature Monitoring 3. Power Control
Maximum Output Power	50 Watts
RF output frequency	450 kHz - 550 kHz
Impedance cut-off	high: 250Ω low: 40Ω

Refer to the RF generator manual for detailed generator operating instructions for RF catheter ablation.

Accessories:

Use only the appropriate Daig extension cable with Redel 14-pin connector to connect to a compatible RF generator.

10 HOW SUPPLIED

The Livewire Catheter is only available with temperature sensing. All catheters use a 4mm tip electrode, Redel 14-pin connector, and standard 2-5-2(mm) spacing. The Livewire TC™ Catheter is available in one of the following six curve styles: Small Curl, Medium Curl, Large Curl, Small Sweep, Medium Sweep & Large Sweep.

10.1 Packaging

The Livewire Catheter is supplied STERILE. The catheter is placed into a thermoformed plastic tray designed to protect the catheter during handling. The tray is placed into a Tyvek® and Mylar® pouch and heat sealed.

10.2 Storage and Shelf-Life

The catheters must be stored in a cool, dark, dry place.
All catheters are labeled with a three year expiration date.

10.3 Warranty and Replacement Policy

Daig Corporation does not warrant or guarantee the Livewire Catheter. Please see the enclosed "Warranty Card" for further information. Also, please see the following labeling sections for information concerning the performance of this device: Contraindications, Warnings and Precautions, and Adverse Events.

DAIG CORPORATION
14901 DeVeau Place
Minnetonka, MN 55345-2126 USA

ST. JUDE MEDICAL EUROPE, INC.
Arianelaan 5
1200 Brussels
Belgium

PH: (800) 328-3873
In MN, call: (612) 933-4700
FAX: (612) 933-0307

PH: 32 2 774 68 11
FAX: 32 2 772 83 84

40841 Rev. Preliminary D
11/19/98

COPYRIGHT 1998

PRINTED IN U.S.A.

DAIG, LivewireTC, the stylized LivewireTC, Fast-Cath and the stylized D
are trademarks of Daig Corporation.

ATAKR® is a registered trademark of Medtronic Inc.

Tyvek® and Mylar® are registered trademarks of E.I. du Pont de Nemours and Company.

U.S. Patent Nos. 5,395,328 & 5,395,329

1 **Information for Patients**
2 **Considering Ablation Treatments**
3 **of Heart Rhythm Problems**



4

5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32

This booklet provides an overview of catheter ablation procedures for the treatment of SVT.

You should consult your physician for a complete explanation of this procedure.

**© 1999
Daig Corporation
14901 DeVeau Place
Minnetonka, MN 55345
800-328-3873 or 612-933-4700**

Daig, Stylized D and Livewire TC are trademarks of Daig Corporation

TABLE of CONTENTS

33
34

35 INTRODUCTION..... 1

36 HOW THE HEART WORKS 1

37 THE PLUMBING..... 2

38 THE WIRING..... 3

39 WHAT HAPPENS WHEN THERE IS SOMETHING WRONG WITH
40 THE WIRING? 5

41 WOLFF-PARKINSON-WHITE (WPW) SYNDROME 6

42 AV NODAL REENTRANT TACHYCARDIA (AVNRT)..... 7

43 RAPID AV NODAL CONDUCTION 7

44 TREATING SVT 8

45 ANTIARRHYTHMIC MEDICATIONS..... 8

46 HEART SURGERY 8

47 CATHETER ABLATION 8

48 WHICH PROCEDURE IS "BEST"? 9

49 HOW DOES CATHETER ABLATION WORK?..... 9

50 WHAT CAN I EXPECT DURING THE PROCEDURE?..... 10

51 WHAT CAN I EXPECT AFTER THE PROCEDURE?..... 14

52 WHAT ARE THE RISKS? 16

53 WHAT ARE THE POTENTIAL BENEFITS? 17

54 GLOSSARY..... 18

55 **INTRODUCTION**

56
57 Your doctor has told you that you may benefit from a
58 catheter ablation procedure. This booklet was designed
59 to provide you with an overview of this procedure and
60 address many of the questions you may have. However,
61 as with all medical procedures, you should consult your
62 physician for a complete explanation of this procedure.

63
64 **HOW THE HEART WORKS**

65
66 Before discussing the specifics of catheter ablation
67 procedures, you will need to know a little about how the
68 heart works.

69
70 The heart is a hollow, muscular organ that pumps blood
71 throughout the body. This is an important function as
72 blood contains both the oxygen and nutrients your body
73 needs to work properly. Blood is pumped from the heart
74 when it contracts (beats). This contraction occurs when
75 the heart muscle is stimulated by an electrical impulse,
76 either naturally or from an external source. To do all of
77 this, your heart has plumbing and wiring, just like your
78 home.
79


80 **The Plumbing**

81

82 Your home no doubt has a plumbing system, which
83 includes sinks, a water heater, pipes, etc. Your heart
84 also has a plumbing system. Instead of sinks, your heart
85 has four compartments, or chambers. These four
86 chambers work together to pump blood throughout the
87 body. The upper two chambers are called the left and
88 right atria (singular: atrium) and the two lower chambers
89 are called the left and right ventricles (see Figure 1).

90

- Figure 1 -



91 The atria receive blood from the body, and the ventricles
92 pump blood out to the body. Normally both atria contract
93 (or pump) at approximately the same time, and then both
94 ventricles contract together a short time later.

95

96 **The Wiring**

97

98 Your home also has an electrical system. So does your
99 heart. The heart depends on its electrical system to
100 stimulate it to contract properly. This system directs
101 electrical impulses through the heart in an organized
102 fashion. Each part of the heart contracts when an
103 impulse reaches it, and each normal heartbeat starts at
104 the sinus node, or SA node (see Figure 2).


105

106 The sinus node is a group of special heart cells located in
107 the right atrium. Normally the sinus node generates an
108 electrical impulse at the proper time for each heartbeat.
109 Therefore, the sinus node is the heart's "natural
110 pacemaker" as it controls the normal heart rate.

111

112 From the sinus node the electrical impulse travels
113 throughout the atria, causing them to contract and
114 squeeze blood into the ventricles.

115

- Figure 2 -

116 From the atria, the electrical impulse travels to the
117 atrioventricular node (or AV node), which is located
118 between the atria and the ventricles. The AV node
119 provides the only normal path for the impulse to travel
120 from the atria to the ventricles. This is because there is a
121 layer of tissue between the atria and ventricles that does
122 not conduct electrical impulses. The AV node sends the
123 impulse to the ventricles at the proper time. This impulse
124 stimulates the ventricles to contract and pump blood
125 efficiently to the body.

126 **WHAT HAPPENS WHEN THERE IS SOMETHING**
127 **WRONG WITH THE WIRING?**


128
129 If the wiring in your house were installed incorrectly it
130 would not function properly. It might even be dangerous.
131 The same is true if there is a problem with the wiring in
132 your heart. Abnormal heart wiring can cause abnormal
133 heart rhythms, or arrhythmia. This usually results in a
134 change in the rate at which the heartbeats.

135
136 An arrhythmia may first be noticed as a skipping or
137 fluttering sensation in the chest (palpitation) or it may
138 cause light-headedness, fainting, chest pain, or
139 shortness of breath. Occasionally, arrhythmias may go
140 unnoticed and are only detected by an electrocardiogram
141 (ECG).

142
143 Some arrhythmias are merely an annoyance. However,
144 others can be serious and can result in heartbeats that
145 are too fast or too slow to pump blood effectively to the
146 body.

147
148 A common type of arrhythmia is called supraventricular
149 tachycardia, or SVT. In this type of arrhythmia the atria
150 beat too fast, often causing the ventricles to also beat too
151 quickly. Three common types or causes of SVT are
152 described in the following pages.

- Figure 3 -



154

155 **Wolff-Parkinson-White (WPW) Syndrome**

156

157 Wolff-Parkinson-White Syndrome is a common cause of
158 SVT. In WPW there is extra tissue connecting the atria
159 and ventricles. This extra tissue is a "short circuit"
160 between these chambers. It provides an extra pathway
161 for electrical impulses to be conducted through the tissue
162 which normally blocks electrical impulses between the
163 atria and ventricles. This short circuit is called an
164 **accessory pathway**, and it allows electrical impulses to
165 travel between the atria and the ventricles without going
166 through the AV node.

167 In WPW an SVT is usually started when an impulse
168 travels down the AV node to the ventricles, and then up
169 through the "short circuit" tissue to the atria. This impulse
170 can then travel through the atria and down the AV node
171 before the sinus node can start the next heartbeat. If the
172 impulse continues to travel in this repeating, circular
173 pattern, it can cause the heart to beat very rapidly.
174

175 **AV Nodal Reentrant Tachycardia (AVNRT)**

176
177 AVNRT is another common form of SVT. In AVNRT
178 there is an extra electrical pathway in or near the AV
179 node. If an electrical impulse is conducted in this
180 pathway, it may direct the impulse through both the AV
181 node and the extra pathway in a repeating, circular
182 pattern. This causes the atria and ventricles to contract
183 with each cycle of this circular pattern, resulting in a rapid
184 heartbeat.
185

186 **Rapid AV Nodal Conduction**

187
188 In some SVTs the atria may spontaneously generate
189 multiple rapid impulses. Many of these impulses can
190 travel through the AV node to the ventricles in an erratic
191 manner. As a result, the heart rhythm can be irregular
192 and rapid. If this happens, the heart will pump blood
193 inefficiently.

194
195
196
197
198
199
200
201
202
203
204
205
206
207
208
209
210
211
212
213
214

TREATING SVT

There are several ways to treat SVT, including antiarrhythmic medications, surgery, and ablation.

Antiarrhythmic medications

change how your heart's electrical system works. They can help prevent abnormal heartbeats from starting or being sustained.

Heart surgery

is rarely performed for the treatment of SVT. During such an operation, surgeons either cut or remove extra electrical pathways causing the arrhythmia.

Catheter ablation

is another procedure in which the extra electrical pathway(s) causing your arrhythmia is disrupted or destroyed.

215 **Which Procedure Is “Best”?**

216

217 You and your doctor should decide which is the best
218 treatment for your particular arrhythmia, after discussing
219 the risks and benefits of each.


220

221 **HOW DOES CATHETER ABLATION WORK?**

222

223 During catheter ablation, doctors insert several special
224 long, flexible tubes with wires called electrode catheters
225 into the heart (see Figure 4). Some of these can be used
226 to study your arrhythmia. However, one of these
227 catheters will be used for the actual ablation. One
228 catheter your doctor may select for this ablation
229 treatment is called the Livewire TC™ ablation catheter
230 and is made by Daig Corporation, Minnetonka,
231 Minnesota. Your physician will position the ablation
232 catheter so that it lies on or very close to the abnormal
233 tissue. High-frequency electrical energy is then delivered
234 through the ablation catheter to this tissue. The small
235 area of heart tissue under the tip of the ablation catheter
236 is heated by this high-frequency energy, destroying the
237 tissue. As a result, this tissue is no longer capable of
238 conducting or sustaining the SVT.

- Figure 4 -



240
241

WHAT CAN I EXPECT DURING THE PROCEDURE?

242 Catheter ablation is done in a specially equipped room,
243 called an Electrophysiology Laboratory (EP Lab).
244 Sometimes these procedures are performed in a Cardiac
245 Catheterization Laboratory (Cath Lab).

246
247 Normally you'll be taken to the EP Lab on a movable bed,
248 then shifted onto a x-ray table. This special table is
249 moveable and has a large x-ray camera over it. There is
250 other equipment in an EP Lab, including viewing screens,
251 heart monitors, emergency equipment, etc. Once on the
252 x-ray table several types of monitoring equipment will be
253 connected to you. You will then be covered with a sterile
254 drape.

255
256 There are usually several EP Lab staff present during the
257 procedure, including one or more electrophysiologists (a
258 cardiologist with special training), nurses, and
259 technicians.

260
261 A nurse will shave and cleanse the area where the
262 catheters will be inserted. In most cases this will be the
263 groin or neck area. To numb the area, a local anesthetic
264 is injected into the skin with a tiny needle.

265
266 A small intravenous needle ("IV line") will be inserted into
267 a vein in your arm. This allows drugs to be injected
268 directly into the vein, if necessary.

269

270 Many times you will be awake during the procedure,
271 although medication is often given to help you relax and
272 be comfortable. However, you may fall asleep during the
273 procedure. The staff will be monitoring you constantly.

274

275 A small incision is made in the numbed skin, again
276 usually in the groin or neck area. A needle is used to
277 puncture the blood vessel (usually a vein, but sometimes
278 an artery) into which an ablation and /or electrode
279 catheter will be inserted.

280

281 One or more electrode catheters are inserted into your
282 blood vessels and advanced toward your heart, while the
283 staff follows the catheter(s) progress on a x-ray viewing
284 screen.

285

286 These electrode catheters can be used to sense
287 electrical activity in various areas of your heart and
288 measure how fast these impulses travel. The electrode
289 catheters can also be used to deliver tiny electrical
290 impulses to stimulate your heart to contract - you should
291 not feel these impulses. By doing so, doctors often
292 attempt to start (or induce) your SVT so they can
293 understand and decide how best to treat it. If you feel the
294 same symptoms you experienced when the arrhythmia

295 occurred previously, you should let the EP Lab staff
296 know.

297

298 Often these induced arrhythmias will stop by themselves.
299 However, if an arrhythmia persists or is very rapid, it may
300 make you faint for a moment. If such an arrhythmia
301 happens, your doctor may decide to deliver an electric
302 shock to your heart in order to stop the abnormal heart
303 rhythm. If you were not in an EP Lab, these arrhythmias
304 could be very dangerous, perhaps even life-threatening.
305 However, well-trained personnel in the EP Lab have the
306 equipment and medications necessary to respond to
307 these arrhythmias.

308

309 The catheter ablation procedure is usually not painful.
310 However, you may feel some pressure at the site(s)
311 where the catheters are inserted. It is also not unusual to
312 experience some mild chest discomfort during the
313 application of the high-frequency energy during the
314 ablation.

315

316 Most catheter ablation procedures are completed within
317 two hours. However, a complete procedure can last up
318 to six hours or more. It is therefore possible that you may
319 feel tired and uncomfortable after lying still for such a
320 lengthy period of time.

321

322 **WHAT CAN I EXPECT AFTER THE PROCEDURE?**

323

324 All the catheters will be removed at the conclusion of the
325 procedure. Firm pressure will be applied at these
326 catheter insertion sites for several minutes in order to
327 prevent bleeding. In some cases a few stitches will be
328 used to close the skin at these locations. A bandage
329 dressing will also be applied.

330

331 Much of the monitoring equipment will be disconnected
332 from you. However, sometimes some of this equipment
333 will remain connected until you have been transported to
334 a recovery area or your hospital room. The IV line in
335 your arm is often left in place.

336

337 Once in the recovery area or your room you will be
338 required to lie flat and still for several hours. You should
339 avoid lifting and bending your leg(s) where the catheters
340 were inserted. This will provide the punctured blood
341 vessels an opportunity to close more completely.

342

343 Typically a nurse will follow your progress for several
344 hours by checking your pulse, blood pressure, and the
345 catheter insertion sites. If you notice bleeding or feel
346 pain at these insertion sites, or if you feel your heart
347 beating rapidly, you should immediately notify the nurse.

348

349 Sometimes you will be allowed to go home the same day
350 as the procedure. However, you may be required to stay
351 in the hospital overnight. Your heart may be monitored
352 with an ECG until you go home. You should make
353 arrangements for someone to take you home from the
354 hospital.

355
356 After you return home you should limit your activity for
357 several days. All vigorous physical exertion and strain
358 (such as lifting heavy objects) should be avoided. In
359 addition, you should carefully follow your doctor's
360 instructions regarding which medications to take.

361
362 You should leave the bandage dressing in place until the
363 next day, or as instructed by your doctor or nurse. They
364 will also tell you how long you should wait after returning
365 home before bathing.

366
367 It is not unusual for there to be a bruise or small lump
368 where the catheters were inserted. This will usually
369 disappear in a week or two. However, it is unusual for
370 this site(s) to become warm to the touch, tender, painful,
371 or for any swelling to increase after you return home. It
372 would also be unusual for you to develop a fever, or
373 experience a recurrence of your rapid heart rhythm, chest
374 pain, dizziness, or shortness of breath. If any of these
375 occur you should contact your doctor immediately.

376

377 **WHAT ARE THE RISKS?**

378

379 Ablation is a procedure that requires the insertion of
380 catheters *into* the body. It therefore does involve some
381 risk.

382

383 Some patients can have bleeding, swelling, or bruising
384 where the catheters were inserted. Serious
385 complications do sometimes occur. These include
386 infection, damage to your heart and/or blood vessels, and
387 blood clots. Deaths are very rare during these
388 procedures.

389

390 It is also possible that the heart's normal electrical
391 system could be damaged during this procedure. If this
392 occurs, an artificial pacemaker implant may be
393 necessary.

394

395 Most patients who undergo catheter ablation do not
396 experience complications. However, you should be
397 aware of the risks. If you should have any questions
398 about the potential risks, you should ask your doctor.

04/27/99

399

400

WHAT ARE THE POTENTIAL BENEFITS?

401

402

Catheter ablation may permanently cure your arrhythmia

403

and may allow you to avoid taking medications while

404

resuming a normal lifestyle.

405

GLOSSARY

406

407

408 Arrhythmia – abnormal heart rhythm.

409

410 Accessory pathway – an abnormal, extra electrical
411 pathway between the upper and lower chambers of
412 the heart.

413

414 Artificial pacemaker - a device often implanted in the
415 body that delivers regular electrical impulses to the
416 heart when the heart's normal electrical system fails.

417

418 Atrium - (pl: atria) upper chamber of the heart. There is
419 one on the right side of the heart and one on the left
420 side. These chambers receive blood from the rest of
421 the body.

422

423 AV node – group of specialized cells that lies between
424 the atria and the ventricles. The heart's electrical
425 impulse normally passes through these cells on its
426 way to the ventricles.

427

428 AV nodal reentrant tachycardia (AVNRT) – very rapid
429 heart rate that is caused by an abnormal pathway
430 close to the AV node.

431

432 Catheter ablation - non-surgical technique that disrupts
433 (destroys) parts of the abnormal electrical pathway
434 that is causing an arrhythmia.

435

436 Electrocardiogram – test that records the electrical
437 activity of the heart. Often called an ECG or EKG.

438

439 Electrophysiology study – heart rhythm study. Often
440 called an EP study.

441

442 Sinus node - group of specialized cells in the right atrium.
443 Normally these cells act as the heart's "natural
444 pacemaker", setting the pace for the heartbeat.

445

446 Supraventricular tachycardia (SVT) – very rapid heart
447 rate that begins in the heart's upper chambers.

448

449 Tachycardia – rapid heart rate.

450

451 Ventricle - lower chamber of the heart. There is one on
452 the right side of the heart and one on the left side.
453 These chambers pump blood out of the heart.

454

455 Wolff-Parkinson-White (WPW) syndrome – a condition in
456 which a very rapid heart rate is caused by an
457 abnormal extra pathway between the upper and lower
458 chambers of the heart.