

**3. 510(K) SUMMARY****Applicant - Manufacturer Name and Address**

Micro Therapeutics dba ev3 Neurovascular  
9775 Toledo Way  
Irvine, CA 92618  
Establishment Registration No. 2029214

FEB - 2 2010

**Date:** January 08, 2010

**Contact Information**

Laura Heaton  
Senior Regulatory Affairs Specialist  
Main: (949) 837-3700  
Direct: (949) 680-1241  
Fax: (949) 859-7228

**Trade Names**

HyperForm™ Occlusion Balloon Catheter  
HyperGlide™ Occlusion Balloon Catheter

**Device Classification & Common Name**

Product Code: MJN, 21CFR Part 870.4450  
Classification: Class II  
Classification Name: Catheter, Intravascular Occluding, Temporary  
Common Name: Occlusion Balloon Catheter

**Predicate Devices - 510(k) References**

HyperGlide™ Occlusion Balloon Catheter	K092495
HyperGlide™ and HyperForm™ Occlusion Balloon Catheters	K091458
HyperGlide™ Occlusion Balloon Catheter	K090728
HyperGlide™ Occlusion Balloon Catheter	K021066
HyperForm™ Occlusion Balloon Catheter	K011656

**Description of the Device Subject to Premarket Notification**

The Occlusion Balloon Catheter is a single lumen balloon catheter that requires the insertion of the 0.010" guidewire to occlude the central lumen to allow inflation of the balloon. When the distal 10 cm platinum coil tip of the guidewire is advanced to or past the catheter tip, it occludes the inflation holes allowing the balloon to inflate through catheter sideholes. The Occlusion Balloon Catheters are marketed as the HyperForm™ and HyperGlide™ Occlusion Balloon Systems which include a guidewire.

**Indications for Use**

The Occlusion Balloon Catheter is designed for the use in the blood vessels of the peripheral and neuro vasculature where temporary occlusion is desired. The Occlusion Balloon Catheter offers a vessel selective technique of temporary vascular occlusion which is useful in selectively stopping or controlling blood flow.

**Performance Data**

The Occlusion Balloon Catheter devices have not changed; therefore there is no performance data included in this submission. The Instructions for Use have been revised to remove the pediatric and neonatal contraindication.

**Substantial Equivalence**

The devices have not changed and no new risks have been identified. The indications for use demonstrate the ev3 Occlusion Balloon Catheters are substantially equivalent to the predicate devices.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room W-066-0609  
Silver Spring, MD 20993-0002

FEB - 2 2010

EV3 Neurovascular  
c/o Laura Heaton  
Senior Regulatory Affairs Specialist  
9775 Toledo Way  
Irvine, CA 92618

Re: K100063

Trade/Device Name:  
Regulation Number: 21 CFR 870.4450  
Regulation Name: Vascular clamp  
Regulatory Class: Class II (two)  
Product Code: MJN  
Dated: January 8, 2010  
Received: January 15, 2010

Dear Ms. Heaton:

We have reviewed your Section 510(k) premarket notification of intent to market the devices referenced above and have determined the devices are substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the devices, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your devices are classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your devices comply with other requirements of the Act


or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your devices on our labeling regulation (21 CFR Part 801), please go to

<http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

K100063

2. STATEMENT OF INDICATIONS FOR USE

*Indications for Use*

510(k) Number (if known):

Device Name: HyperForm™ Occlusion Balloon Catheters

Indications for Use: The HyperForm™ Occlusion Balloon Catheters are indicated for the use in blood vessels of the peripheral and neuro vasculature where temporary occlusion is desired. These catheters offer (1) a vessel selective technique of temporary vascular occlusion, which is useful in selectively stopping or controlling blood flow, and (2) for balloon-assisted embolization of intracranial aneurysms.

Prescription Use  X

AND/OR

Over-The-Counter Use \_\_\_\_\_

(Part 21 CFR 801 Subpart D)

(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of 1

  
\_\_\_\_\_  
(Division Sign-Off)  
Division of Cardiovascular Devices  
510(k) Number K100063