

K100403

510(k) SUMMARY

JUN 28 2010

Submitted by: Masimo Corporation
 40 Parker
 Irvine, CA 92618
 949-297-7000
 FAX 949-297-7592

Company Contact: Marguerite Thomlinson, Manager of Regulatory Affairs

Date Summary Prepared: April 23, 2010

Trade Name Masimo Rainbow SET Pronto-7 Pulse CO-Oximeter and Accessories

Common Name Oximeter

Classification Name/ Product Code: Oximeter Section 870.2700/ Product Code DQA

Predicate Device: Masimo Rainbow SET[®] Pronto Pulse CO-Oximeter and Accessories, 510(k) Number – K091057

Device Description

The Masimo Rainbow SET Pronto-7 Pulse CO-Oximeter and Accessories (Pronto-7) include the noninvasive technology of measurement for functional saturation of arterial oxygen hemoglobin (SpO₂), pulse rate and total hemoglobin concentration (SpHb), which is based on the same technology of the predicate device, the Masimo Rainbow SET[®] Pronto Pulse CO-Oximeter and Accessories (Pronto), 510(k) no. K091057.

Intended Use/ Indications for Use

The Masimo Rainbow SET Pronto-7 Pulse CO-Oximeter and Accessories are indicated for noninvasive spot checking of functional saturation of arterial oxygen hemoglobin (SpO₂), pulse rate, and total hemoglobin concentration (SpHb). The Masimo Rainbow SET Pronto-7 Pulse CO-Oximeter and Accessories are indicated for use by trained personnel, with adult and pediatric individuals, in clinical and non-clinical settings (e.g., hospitals, hospital-type facilities, home, clinics, physician offices, blood donation facilities, and ambulatory surgery centers).

Technology Comparison

Features in the Pronto-7, which are not in the Pronto, include touch-screen user interface, pleth waveform display, rechargeable batteries or AC power, USB connection, flash memory, earphone jack, and wireless connection. The Pronto-7 has the following specifications:

FEATURES	SPECIFICATIONS
Display Range	Saturation (SpO ₂): 0-100%; Pulse Rate (bpm): 30-250 bpm Total Hemoglobin (SpHb): 0-25 g/dl; Perfusion Index (PI): 0.02-20%
Accuracy: Adults and Pediatrics > 30kg	SpO ₂ : 70-100±2%; Pulse Rate: 30-250±3 bpm; SpHb: 6-18 g/dl ±1 g/dl
Resolution	SpO ₂ : 1%; Pulse Rate: 1 bpm; SpHb: 0.1 g/dl
AC Power	Voltage Input: 100-240 VAC, 50-60 Hz, Max Power Consumption: 15 VA
Batteries	Rechargeable lithium polymer

510(k) SUMMARY

FEATURES	SPECIFICATIONS
Temperature	Operating/ Storage: 41 to 104°F (5 to 40°C)/ -40 to 158°F (-40 to 70°C)
Relative Humidity	5 to 95% noncondensing
Operating Altitude	500 to 1,060 mbar pressure; -1,000 to 18,000 ft (-304 to 5,486m)
Alarm	System Failure; Low Battery Alarm
Display and Indicators	SpO ₂ (%); Pulse Rate (bpm); SpHb (g/dl); Perfusion Index (%) Pleth Waveform; Sensor Status; Status Messages; Battery Status
Connection/ Output	Wireless: Bluetooth 2.0; WiFi b/g; Flash Memory: Micro SD card slot
Earphone jack	Connection for standard 3.5 mm earphone jack
EMC/ Electrical Safety Compliance	EN 60601-1-2, Class B IEC 60601-1, UL 60601-1, Internally Powered, AC Power Class 2 Type BF-Applied Part, IPX1, Class 2
Mode of Operation	Spot check

Clinical Data Summary

Clinical Studies: SpO₂ accuracy has been validated in human blood studies on healthy adult male and female volunteers with light to dark skin pigmentation. Pronto-7 measurements on volunteers induced hypoxia in the range of 70-100% SpO₂ are compared against a laboratory CO-oximeter and ECG monitor. SpHb accuracy has been validated with (arterial/ venous) blood from healthy adult male and female volunteers and on patients with light to dark skin pigmentation in the range of 6-18 g/dl SpHb, with Pronto-7 measurements compared against a laboratory CO-oximeter. Pulse rate accuracy has been validated on healthy adult male and female volunteers and on patients with light to dark skin pigmentation in the range of 40-110 bpm. The variation in these studies equals plus or minus one standard deviation encompasses 68% of the population.

Clinical Results: No device-related adverse events.

The clinical studies were performed in accordance with ISO 9919:2005. The studies resulted in SpO₂ accuracy (rms) ≤ 2%, pulse rate accuracy (rms) ≤ 3bpm, and SpHb accuracy (rms) ≤ 1g/dl.

Non-Clinical Data Summary

The Pronto 7 complies with the voluntary standards as detailed in this submission. Laboratory testing for biocompatibility, safety and environmental was conducted to verify that the Pronto 7 met all design specifications and was substantially equivalent to the predicate device.

Conclusions

The information in this 510(k) submission demonstrates that the Masimo Rainbow SET Pronto 7 Pulse CO-Oximeter and Accessories are substantially equivalent to the predicate device, with respect to safety, effectiveness, and performance.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room - WO66-G609
Silver Spring, MD 20993-0002

Ms. Marguerite Thomlinson
Manager of Regulatory Affairs
Masimo Corporation
40 Parker
Irvine, California 92618

JUN 23 2010

Re: K100403
Trade/Device Name: Masimo Rainbow SET Pronto-7 Pulse CO-Oximeter and
Accessories
Regulation Number: 21 CFR 870.2700
Regulation Name: Oximeter
Regulatory Class: II
Product Code: DQA
Dated: June 11, 2010
Received: June 14, 2010

Dear Ms. Thomlinson:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to

<http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Anthony D. Watson, BS, MS, MBA
Director
Division of Anesthesiology, General Hospital
Infection Control and Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): _____

Device Name: Masimo Rainbow SET Pronto-7 Pulse CO-Oximeter and Accessories

Indications For Use:

The Masimo Rainbow SET Pronto-7 Pulse CO-Oximeter and Accessories are indicated for noninvasive spot checking of functional saturation of arterial oxygen hemoglobin (SpO₂), pulse rate, and total hemoglobin concentration (SpHb). The Masimo Rainbow SET Pronto-7 Pulse CO-Oximeter and Accessories are indicated for use by trained personnel, with adult and pediatric individuals, in clinical and non-clinical settings (e.g., hospitals, hospital-type facilities, home, clinics, physician offices, blood donation facilities, and ambulatory surgery centers).

Prescription Use X
(Per 21 CFR 801.109 Subpart D)

AND/OR

Over-The-Counter Use _____
(Per 21 CFR 801.109 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF
NEEDED)



Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)
Division of Anesthesiology, General Hospital
Infection Control, Dental Devices

510(k) Number: K100403