

K100429

1/2

5. 510(k) Summary

Additional Information to Premarket Notification Summary

JUN 11 2010

1. Sponsor Information:
Laboratories URGO
42 Rue de Longvic
21300 Chenove
France

Contact Person: Sophie Fortin
Regulatory Affairs Director
Phone: +33.(0)3.80.44.28.78
Fax: +33.(0)3.80.44.71.40
Email: s.fortin@fr.urgo.com

2. Device Name:
Common or Usual Name: Antimicrobial Wound Dressing with Silver
Proprietary Name: Urgocell® Ag
Classification Name: Unclassified
3. Predicate Devices:
Urgotul Ag (K061220), Laboratoires URGO
Urgocell Ag (K062559), Laboratoires URGO
4. Description of Device
Urgocell® Ag, antimicrobial wound dressing with silver, is a non-adhesive, non-occlusive, antimicrobial absorbent dressing, composed of 3 layers:
 - in contact with the wound, a polyester mesh impregnated with a matrix comprising of hydrocolloid particles (carboxymethylcellulose), cohesion polymers, petroleum jelly and silver sulphate (3.22 mg/sq.in)
 - a non-sensitising, super-absorbent polyurethane foam pad,
 - a protective, semi-permeable polyurethane backing.
5. Indications for Use
The barrier functions of Urgocell® Ag, antimicrobial wound dressing with silver, may help reduce infection in moderately to high exuding partial and full thickness wounds, including decubitus ulcers, venous static ulcers, diabetic ulcer, first and second degree burns, donor and graft sites.

2/2

6. Technological Characteristics

FUNCTION	Urgocell Ag (K062559)	Urgocell Ag
Indications for Use	The barrier functions of Urgocell® Ag, antimicrobial wound dressing with silver, may help reduce infection in moderately to high exuding partial and full thickness wounds, including decubitus ulcers, venous static ulcers, diabetic ulcer, first and second degree burns, donor and graft sites.	The barrier functions of Urgocell® Ag, antimicrobial wound dressing with silver, may help reduce infection in moderately to high exuding partial and full thickness wounds, including decubitus ulcers, venous static ulcers, diabetic ulcer, first and second degree burns, donor and graft sites.
Structure	Composed of three layers: - in contact with the wound, a polyester mesh impregnated with a matrix comprising of hydrocolloid particles (carboxymethylcellulose), cohesion polymers, petroleum jelly and silver sulphate, - a non-sensitising, super-absorbent polyurethane foam pad, - a protective, semi-permeable polyurethane backing	Composed of three layers: -in contact with the wound, a polyester mesh impregnated with a matrix of hydrocolloid particles (carboxymethylcellulose), cohesion polymers, petroleum jelly and silver sulphate, - a non-sensitising, super-absorbent polyurethane foam pad, - a protective, semi-permeable polyurethane backing
Antimicrobial activity	Studies of antimicrobial activity against MRSA, Pseudomonas aeruginosa, and Streptococcus pyogenes	Studies of antimicrobial activity against MRSA, Pseudomonas aeruginosa, Streptococcus pyogenes, Enterococcus faecalis, Escherichia coli, and Candida albicans

7. Description of Safety and Substantial Equivalence:

Safety Studies

The standard battery of safety and biocompatibility studies that were previously conducted showed the Urgocell Ag to be comparable to predicate devices and other silver containing wound dressings. These studies included: Cytotoxicity, Irritation study in rabbit, and Sensitization study in guinea pigs.

Studies of the antimicrobial activity of Urgocell Ag wound dressing (K062559) were previously conducted against MRSA, Pseudomonas aeruginosa, and Streptococcus pyogenes. With this submission, data was submitted concerning the antimicrobial activity of Urgocell® Ag on three new strains: *Enterococcus faecalis*, *Escherichia coli* and *Candida albicans*.

Substantial Equivalence

Urgocell Ag compared to Urgocell Ag (K062559) shows substantial equivalence in comparable functions and indications for use.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room -WO66-G609
Silver Spring, MD 20993-0002

JUN 11 2010

Laboratories URGO
% Ms. Sophie Fortin
Regulatory Affairs Director
42 Rue de Longvic
21300 Chenove, France

Re: K100429

Trade/Device Name: Urgocell[®] Ag Antimicrobial Wound Dressing with Silver
Regulatory Class: Unclassified
Product Code: FRO
Dated: May 25, 2010
Received: June 03, 2010

Dear Ms. Fortin:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set


Page 2 - Ms. Sophie Fortin

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Mark N. Melkerson
Director
Division of Surgical, Orthopedic
And Restorative Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

4. Statement of Indications for Use

510(k) Number (if known): K100429

Device Name: URGOCELL AG, Antimicrobial Wound Dressing with Silver

Indications For Use:


The barrier functions of Urgocell Ag, antimicrobial wound dressing with silver, may help reduce infection in moderately to high exuding partial and full thickness wounds, including decubitus ulcers, venous static ulcers, diabetic ulcer, first and second degree burns, donor and graft sites.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR Over-The-Counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)
Division of Surgical, Orthopedic,
and Restorative Devices

Page 1 of 1

510(k) Number K100429