K100447

Premarket Notification

Section 5: 510(k) Summary

Device Information:			
Category	Comments		
Sponsor / Submitter:	MyoScience, Inc		
	525 Chesapeake Drive		
	Redwood City, CA 94063		
	(650) 474-2600		
	(650) 474-2700		
Correspondent Contact	Tracey Henry		
Information:	Sr. Director RAQA, Clinical Compliance		
	525 Chesapeake Drive		
	Redwood City, CA 94063		
	(650) 474-2600		
	(650) 474-2900		
Device Common Name:	Cryogenic Surgical device		
Device Classification & Code:	Class II, GXH		
Device Classification Name:	Cryosurgical unit and accessories (21 CFR 882.4250)		
Device Proprietary Name:	MyoScience Cryo-Touch		

a. Predicate Device Information:

The Cryo-Touch is substantially equivalent to the following currently legally marked devices:

510(k) Number	Product	Sponsor
K083493	Cryo-Touch	MyoScience, Inc
K050272	Cryo-PaC [™] System	Cryomedical Instruments
K854334	Painblocker WA 5000	Wallach Surgical Devices
K831963	Cryomedics Neurostat	Spembly Medical Ltd

b. Date Summary Prepared

June 21, 2010

c. Description of Device

The MyoScience Cryo-Touch is a hand-held, single patient-use disposable, cryogenic device used to destroy nervous tissue or produce lesions in peripheral nervous tissue by the application of cold to the selected site for the purposes of blocking pain. The Cryo-Touch has been developed for applications in cryoanalgesia. Cryoanalgesia for intractable pain involves the location and freezing of the nerve associated with the pain.

The device is based on introduction of a cryoprobe cooled by the cryogenic fluid (liquid nitrous oxide (N_2O)) to the selected site and activation of the freeze control. An iceball forms around the tip of the needle assembly and the adjacent nerve.

The Cryo-Touch may be used with a standard off-the-shelf nerve stimulator device in applications where precise nerve location is desired.

Device Design

The device is comprised of three main components:

Page 1

Section 5: 510k Summary

K100447

Premarket Notification

- 1) a control unit (controller) with cryogen cap,
- 2) needle assembly (Cryoprobe), and
- 3) a cryogen cylinder(Nitrous Oxide).

The entire device is hand-held and for single patient-use.

The Cryo-Touch controller is battery operated; the battery powers electronics within the controller that control the nitrous oxide flow and regulate light emitting diodes (LEDs) that provide feedback to the user when the device is ready to use.

An assortment of needle assemblies is available for the Cryo-Touch, ranging from a single to a dual needle, from 25 gauge to 30 gauge, and from 10mm to 25 mm in length. All needle assemblies are made of stainless steel and have a closed-tip that fully contains the cryogen so that it does not enter the target tissue. The needle assembly is the only patient contacting component of the Cryo-Touch. The user removes the needle assembly from the sterile packaging and attaches it to the controller by inserting the needle hub end into the distal end of the controller.

The cryogen is provided as a nitrous oxide cylinder. The cartridge is filled with liquid N_2O (83%) and the rest with N_2O gas. A single nitrous oxide cylinder treats for approximately 180 seconds. The cylinder, once empty, may be replaced with a new cylinder by the user.

Device Functionality/Scientific Concepts

The device functionality is based on the user introducing the needle assembly proximate to the target nervous tissue. The user then initiates the flow of cryogen by opening the start/stop valve. Liquid cryogen flows from the controller into the closed-tip needle assembly. The needle is cooled by the Joule-Thompson effect; as the liquid cryogen expands into a gas, an iceball develops around the external surface of the needle assembly causing the surrounding tissue to be frozen. As the iceball grows, its leading edge advances through the tissue and tissue that comes into contact with the needle is frozen to a nominal target temperature of -65°C (-85°F). The user then stops the flow of cryogen by closing the start/stop valve and allows the iceball to thaw before removing the needle assembly.

d. Intended Use

The MyoScience Cryo-Touch is used to produce lesions in peripheral nervous tissue by the application of cold to the selected site for the blocking of pain. The Cryo-Touch is <u>not</u> indicated for treatment of central nervous system tissue.

e. Comparison to Predicate Devices

The MyoScience Cryo-Touch is substantially equivalent in intended use, technology, design and materials to the above listed legally marketed predicate devices. The technological characteristics of the new device in comparison to the predicate devices are summarized below:

Page



••••

Premarket Notification

·····	This Device	Predicate Device				
Parameter	MyoScience	MyoScience	Cryomedical	Wallach Surgical	Spembly	
	Cryo-Touch	Cryo-Touch	Instruments	Devices	Medical Ltd	
		ко83493	Cryo-PaC™	Painblocker	Cryomedics	
			System K050272	WA5000 K854334	Neurostat	
					K831963	
Intended Use	Identical to all predicates	Destroy tissue through freezing				
Target	Identical to all	Aduits				
Population	predicates					
Anatomical Sites	Identical to predicates	Not specified	pecified Peripheral nerves			
Intended	identical to all	Qualified medical personnel (doctors, specialists) in hospital or medical environment				
users/clinical	predicates					
setting						
Technology	Identical to Cryo-	Cryogenic surgical device with needle which penetrates treatment area.				
	Touch. Similar to					
	Cryo-Pac, Baiablacker WAS000					
	and Neurostat.					
Energy used/or	Identical to all	Cryotherapy/	Cryotherapy/	Cryotherapy/	Cryotherapy/	
delivered	predicâtes	Nitrous Oxide	Nitrous Oxide	Nitrous Oxide or	Nitrous Oxide or	
	· · · · · · · · · · · · · · · · · · ·		· · · · · · · · · · · · · · · · · · ·	Carbon Dioxide	Carbon Dioxide	
Human Factors	Identical to Cryo-	Hand-held and or portable device containing cryogen. Some units have footswitch				
	Touch.	and detachable cryoprobes.				
Biocompatibility	Identical to all predicates	Biocompatible patient contacting materials				
Operating	Identical to all	Joule-Thomson Effect				
Principle	predicates					
Patient	Identical to Cryo-	Closed-tip stainless steel cryoprobe				
contacting	Touch	Ranging from 12 – 30 gauge				
materials						
Treatment	Identical to Cryo-Pac,	-15° C to -25°C				
Temperature	Painblocker WA5000	(5°F to -13°F) for Cryo-Touch, -55°C to -80 (-67°F to -112°F) for other predicates				
	and Neurostat					
Power Source	Yes. Identical to	Battery powered	Mains powered	8attery powered	Battery powered	
1	Cryo-Touch			l		

The device that is the subject of this 510(k) application is identical to the previously cleared Cryo-Touch device with the following exceptions:

- (1) A proprietary filter was added internally. The filter has no direct or indirect patient contact.
- (2) The Cryo-Touch now delivers a target treatment temperature of -65°C which is within the range of the other predicate devices, the Cryo-Pac, the Painblocker WA5000 and the Neurostat.

f. Summary of Supporting Data

Nonclinical testing Submitted: Verification testing was performed on the Cryo-Touch device to demonstrate that the product met the design requirements for system performance and temperature profiling during simulated use conditions.

Page <u>3</u>

K100447

Premarket Notification

A preclinical animal study designed to evaluate the effects of cryotreatment with the Cryo-Touch device on peripheral nerves in a rat model was performed. The study evaluated the impact of the cryolesioning treatment on peripheral nerves in terms of physiologic function and histologic changes.

All other relevant performance testing was submitted as part of the submission for the predicate Cryo-Touch device.

Clinical Testing Submitted: None

g. Conclusion

MyoScience concludes that the Cryo-Touch described in this submission is substantially equivalent to the predicate devices.



Food and Drug Administration 10903 New Hampshire Avenue Document Control Room –WO66-G609 Silver Spring, MD 20993-0002

MyoScience, Inc. c/o Ms. Tracey Henry Sr. Director, Regulatory Affairs and Quality Assurance 525 Chesapeake Drive Redwood City, CA 94063

JUN 2 5 2010

Re: K100447

Trade/Device Name: Cryogenic Surgical Device Regulation Number: 21 CFR 882.4250 Regulation Name: Cryosurgical unit and accessories Regulatory Class: II Product Code: GXH Dated: June 21, 2010 Received: June 22, 2010

Dear Ms. Henry:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

Page 2 - Ms. Tracey Henry

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <u>http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm</u> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm.

Sincerely yours,

esia Alexander

Malvina B. Eydelman, M.D.
Director
Division of Ophthalmic, Neurological and Ear, Nose and Throat Devices
Office of Device Evaluation Health
Center for Devices and Radiological Health

Enclosure

*MyoScience

Cryo-Touch Premarket Notification

Section 4: Indications for Use Statement

510(k) Number:

K100447

Device Name:

MyoScience Cryo-Touch

Indications for Use:

The MyoScience Cryo-Touch is used to produce lesions in peripheral nervous tissue by the application of cold to the selected site for the blocking of pain. The Cryo-Touch is <u>not</u> indicated for treatment of central nervous system tissue.

Prescription Use <u>X</u> (Part 21 CFR 801 Subpart D) AND/OR Over-the-Counter Use _____ (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of _____

Kautma

(Division Sign-Off) Division of Ophthalmic, Neurological and Ear, Nose and Throat Devices

1*004*47 510(k) Number

Page <u>1</u>

Section 4: Indications for Use Statement