

K102008

510(k) Summary

FDA CDRH DMC

A. Device Name

Proprietary Name	Glidesheath
Classification Name	Catheter Introducer (as per 870.1340)
Common Name	Introducer Sheath
Product Code	DYB

JUL 21 2010

Received

B. Intended Use

The Glidesheath is used to facilitate placing a catheter through the skin into a vein or artery including but not limited to the radial artery.

The Entry Needle is an accessory device which is used to gain access to the vein or artery including but not limited to the radial artery, for placement of the Mini Guide Wire.

The Mini Guide Wire is an accessory device which is used for placement of the sheath into the vein or artery including but not limited to the radial artery.

C. Device Description

The Glidesheath is comprised of an introducer sheath and a dilator. The Glidesheath is coated with a hydrophilic coating to reduce the frictional resistance of the sheath when inserting or removing the sheath from the patient's blood vessel. The Sheath and Dilator contain bismuth, making these devices visible under fluoroscopy. The Glidesheath is used to facilitate placing a catheter through the skin into a vein or artery including but not limited to the radial artery.

The Entry Needle is an accessory device which is used to gain access to the vein or artery including but not limited to the radial artery, for placement of the Mini Guide Wire.

The Mini Guide Wire is an accessory device which is used for placement of the sheath into the vein or artery including but not limited to the radial artery. The Guide Inserter which is attached to the Mini Guide Wire holder is used to straighten out the wire.

Accessories to the Glidesheath are the metal entry needle and the mini guide wire. Both the metal entry needle and the mini guide wire are packaged with the Glidesheath in a pouch prior to sterilization.

D. Principle of Operation / Technology

The Glidesheath and its accessories are operated manually or by a manual process.

Using the Guide Inserter to straighten the Mini Guide Wire, the Mini Guide Wire is inserted through a cannula placed in the patient's blood vessel. The Glidesheath is then inserted over the Mini Guide Wire and into the blood vessel. The Mini Guide Wire is then withdrawn from the vessel. The Dilator maintains the integrity of the Sheath and dilates the blood vessel while the Glidesheath is being placed into the vessel. The Dilator can be removed and an appropriate catheter can then be inserted.

E. Design / Materials

	Current Glidesheath	Predicate Glidesheath (K082644)
Sheath		
Tubing	Ethylene-Tetrafluoroethylene (ETFE) copolymer containing bismuth trioxide	Ethylene-Tetrafluoroethylene (ETFE) copolymer containing bismuth trioxide
Hydrophilic Coating	Dimethyl acrylamide-glycidyl methacrylate copolymer	Dimethyl acrylamide-glycidyl methacrylate copolymer
Housing	Polypropylene	Polypropylene
Valve	Silicone Rubber	Silicone Rubber
Cap	Polypropylene	Polypropylene
Caulking Pin	Stainless Steel	Stainless Steel
Sheath Support	Styrene elastomer	Styrene elastomer
Side Tube	Polybutadiene	Polybutadiene
Three-way stopcock		
Holder	Polycarbonate	Polycarbonate
Cock	Polyethylene	Polyethylene
Luer Cap	Polypropylene	Polypropylene
Dilator		
Tube	Polypropylene containing bismuth subcarbonate	Polypropylene containing bismuth subcarbonate
Hub	Polypropylene	Polypropylene

	Current Glidesheath	Predicate Glidesheath (K082644)	
Caulking Pin	Stainless Steel	Stainless Steel	
Mini Guide Wire		Plastic Option	Metallic Option
Spring	Palladium	Polyurethane/tungsten with silicon coating	Stainless Steel
Core	Nickel-Titanium alloy	Nickel-Titanium alloy	Stainless Steel
Insertor	Polyethylene	Polyethylene	
Entry Needle			
Needle Hub	Polycarbonate	Polycarbonate	
Cannula	Stainless Steel	Stainless Steel	
Cap	Polypropylene	Polycarbonate	

Both devices have similar parts which function in the same manner. Differences in materials between the Glidesheath device covered in this submission and the predicate device, cleared under K082644, raise no new issues of safety or effectiveness.

F. Specifications

	Current Glidesheath	Predicate Glidesheath (K082644)
Sheath Size	5 & 6 French	4, 5 & 6 French
Sheath Length	10 cm	10 – 25 cm
Sheath Hydrophilic Coating	10 cm (entire length of shaft)	10 – 25cm (entire length of shaft)
Dilator Length	15.5 cm	15.5 – 30.5 cm
Guide Wire OD	0.021”	0.021” – 0.038”
Guide Wire Length	45 cm	10 – 180 cm
Entry Needle: Size	21G	20 G – 21 G
Entry Needle Length	1.5”	1.5”
IV Catheter	Not Available	16 G – 22 G
IV Catheter Length	Not Available	1” – 2.5”
Scalpel	Not Available	Available
Syringe	Not Available	Available
Obturator	Not Available	Available

All specifications of the Glidesheath included in this submission are within the specification ranges of the predicate device. The comparison of specifications raises no new issues of safety or effectiveness.

G. Performance

The Glidesheath successfully passed all of the following performance tests:

Needle

- Needle surface free from defects
- Needle OD
- Needle length
- Needle ID
- Needle hub conical entry angle
- Bevel indicator visibility
- Bevel indicator position
- Needle to hub joint strength
- Gauge Luer taper
- Liquid leakage from fitting assembly under pressure
- Air leakage into the fitting assembly during aspiration
- Separation force of fitting assembly
- Unscrewing torque of fitting assembly
- Ease of assembly
- Resistance to overriding
- Stress cracking
- Corrosion resistance

Guide Wire

- Guidewire surface free from defects
- Tip buckling test
- Test for resistance of guidewires to damage by flexing
- Test for fracture of guidewires
- Test for distal tip retention
- Guidewire OD
- Guidewire length
- Test for corrosion resistance

Dilator

- Dilator surface free from defects
- Dilator tip ID
- Dilator to hub joint strength

Sheath

- Sheath surface free from defects
- Sheath tip ID
- Sheath to housing joint strength

- Housing to cap joint strength

Inserters

- Guidewire inserter surface free from defects

System

- System use in model

H. Biocompatibility and Sterilization

The Glidesheath is classified as Externally Communicating Devices, Circulating Blood, Limited Contact (≤ 24h). Results of the testing demonstrate that the blood contacting materials are biocompatible.

Blood contacting materials were tested in accordance with the tests recommended in the FDA General Program Memorandum #G95-1 (5/1/95): Use of International Standard ISO-10993, “Biological Evaluation of Medical Devices Part-1: Evaluation and Testing”. The Glidesheath successfully passed all of the following biocompatibility tests:

- Physicochemical Profile
- Cytotoxicity
- Sensitization
- Acute Intracutaneous Reactivity
- Acute Systemic Toxicity
- Hemolysis
- Thrombogenicity
- Complement Activation Assay
- Unactivated Partial Thromboplastin Time Assay
- In Vitro Hemolysis
- Genotoxicity
- Pyrogen Study
- Extractable Metals and Acidity/Alkalinity (needle assembly)

Sterilization conditions have been validated according to ANSI / AAMI / ISO 11135, *Sterilization of Health Care Products– Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices* to provide a Sterility Assurance Level (SAL) of 10⁻⁶.

I. Substantial Equivalence

The performance of the Glidesheath in this submission is substantially equivalent to the performance of the predicate device. The equivalence was shown through

comparison of component materials and specifications, performance and biocompatibility testing and sterilization validation.

The Glidesheath is substantially equivalent in intended use, design, technology/principles of operation, materials and performance to the predicate device the Glidesheath, cleared under K082644. Differences between the devices do not raise any significant issues of safety or effectiveness.

J. Submitter Information

Prepared By: Mr. Daniel R. Plonski
Regulatory Affairs Specialist

Prepared For: Terumo Medical Corporation
950 Elkton Blvd.
Elkton, MD 21921
Phone: (410) 392-7213
Fax: (410) 398-7395
Email: daniel.plonski@terumomedical.com

Date Prepared: July 2, 2010



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room - WO66-G609
Silver Spring, MD 20993-0002

Terumo Medical Corporation
c/o Mark Job
Reviewer
Regulatory Technology Services LLC
1394 25th Street NW
Buffalo, MN 55313

JUL 21 2010

Re: K102008
Trade/Device Name: Glidesheath
Regulation Number: 21 CFR 870.1250
Regulation Name: Catheter, Introducer
Regulatory Class: Class II (Two)
Product Code: DYB
Dated: July 15, 2010
Received: July 16, 2010

Dear Mr. Job:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Page 2 – Mr. Mark Job

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



 Bram D. Zuckerman, M.D.
Director
Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

K102008

Indications for Use

510(k) Number (if known): K102008

Device Name: Glidesheath™

Indications For Use:

JUL 21 2010

The Glidesheath is used to facilitate placing a catheter through the skin into a vein or artery including but not limited to the radial artery.

The Entry Needle is an accessory device which is used to gain access to the vein or artery including but not limited to the radial artery, for placement of the Mini Guide Wire.

The Mini Guide Wire is an accessory device which is used for placement of the sheath into the vein or artery including but not limited to the radial artery.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Diana R. Vidmer

Division Sign-Off)
Division of Cardiovascular Devices

510(k) Number K102008