



GE Healthcare
510(k) Premarket Notification Submission

510(k) Summary

OCT 20 2010

In accordance with 21 CFR 807.92 the following summary of information is provided:

Date: August 12, 2010

Submitter: GE Medical Systems *Information Technologies, Inc.*
8200 Tower Avenue
Milwaukee, WI 53223

Primary Contact Person: Mary Carter
Regulatory Affairs Leader
GE Healthcare
414-362-2626
414-362 2585

Secondary Contact Person: David Wahlig
Director, Regulatory Affairs
GE Healthcare
414-362-3242
414-362 2585

Device: Trade Name: CARESCAPE V100 Vital Signs Monitor

Common/Usual Name: Physiological or Vital Signs Monitor, Patient Monitor

Classification Names: 21 CFR 870.2300

Product Code: MWI

Predicate Device(s): K073203 CARESCAPE V100
K011291 Exergen TemporalScanner Thermometer, SensorTouch
K050230 TurboTemp Trio

Device Description: The CARESCAPE V100 vital signs monitor provides a small, portable, easy-to-use monitoring alternative for sub-acute hospital and non-hospital settings. The monitor is for use on adult, pediatric, or neonatal patients—one at a time. The battery-operated monitor offers noninvasive determination of systolic blood pressure, diastolic blood pressure, mean arterial pressure, pulse rate, oxygen saturation, and temperature. Monitors are available with or without integrated printers as well as the following parameters and technologies.

- NIBP, Pulse: SuperSTAT, Auscultatory, or Classic
- SpO2: Ohmeda TruSignal, Nellcor OxiMax, or Masimo SET
- Temperature: Alaris Turbo Temp, Alaris Tri-Site, or Exergen TemporalScanner Thermometer

Intended Use: The CARESCAPE V100 Vital Signs Monitor is intended to monitor a single adult, pediatric or neonatal patient's vital signs at the bedside or during intra-hospital transport. Vital signs parameters include non-invasive blood pressure (systolic, diastolic, and mean arterial pressure), pulse rate, and/or oxygen saturation (pulse oximetry) and/or temperature. The portable device is designed for use in numerous clinical settings in various hospital departments such as emergency, radiology, recovery, medical/surgical, labor and delivery, endoscopy, cardiac step-down. The CARESCAPE V100 Vital Signs Monitor can also be used in satellite areas, physicians' office, or alternate care settings.



GE Healthcare
510(k) Premarket Notification Submission

Technology: The CARESCAPE V100 Vital Signs Monitor employs the same fundamental scientific technology as its predicate devices.

Determination of Substantial Equivalence: Summary of Non-Clinical Tests:
The CARESCAPE V100 Vital Signs Monitor and its applications comply with voluntary standards as detailed in Section 9, 11 and 17 of this premarket submission. The following quality assurance measures were applied to the development of the system:

- Risk Analysis
- Requirements Reviews
- Design Reviews
- Testing on unit level (Module verification)
- Integration testing (System verification)
- Performance testing (Verification)
- Safety testing (Verification)
- Simulated use testing (Validation)

Summary of Clinical Tests:

The subject of this premarket submission, CARESCAPE V100 Vital Signs Monitor, did not require clinical studies to support substantial equivalence.

Conclusion: GE Healthcare considers the CARESCAPE V100 Vital Signs Monitor to be as safe, as effective, and performance is substantially equivalent to the predicate device(s).



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room –WO66-G609
Silver Spring, MD 20993-0002

Ms. Mary Carter
Regulatory Affairs Leader
GE Healthcare
GE Medical Systems Information Technologies, Inc.
8200 Tower Avenue
Milwaukee, WI 53223

OCT 20 2010

Re: K102426.
Device Name: CARESCAPE V100 Vital Signs Monitor
Regulation Number: 21 CFR 870.2300
Regulation Name: Patient physiological monitor (without arrhythmia detection or alarms)
Regulatory Class: Class II (Two)
Product Codes: MWI, DQA, DXN, FLL
Dated: September 20, 2010
Received: October 8, 2010

Dear Ms. Carter:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Bram D. Zuckerman, M.D.

Director

Division of Cardiovascular Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure



GE Healthcare

510(k) Premarket Notification Submission

510(k) Number (if known): K102426

OCT 20 2010

Device Name: CARESCAPE V100 Vital Signs Monitor

Indications for Use:

The CARESCAPE V100 Vital Signs Monitor is intended to monitor a single adult, pediatric or neonatal patient's vital signs at the bedside or during intra-hospital transport. Vital signs parameters include non-invasive blood pressure (systolic, diastolic, and mean arterial pressure), pulse rate, and/or oxygen saturation (pulse oximetry) and/or temperature.

The portable device is designed for use in numerous clinical settings in various hospital departments such as emergency, radiology, recovery, medical/surgical, labor and delivery, endoscopy, cardiac step-down. The CARESCAPE V100 Vital Signs Monitor can also be used in satellite areas, physicians' office, or alternate care settings.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use
(Part 21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)
Division of Cardiovascular Devices

510(k) Number K102426

1071