

## 510(k) Summary

JAN 20 2011

K102814

---

**Introduction** According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.

---

**Submitter Name, Address, Contact** Roche Diagnostics  
9115 Hague Road  
Indianapolis, IN 46250

Contact person: Kelly French  
Phone: (317) 521-3208  
Fax: (317) 521-2324  
Email: [kelly.french@roche.com](mailto:kelly.french@roche.com)

Secondary contact: Stephanie Greeman  
Phone: (317) 521-2458  
Fax: (317) 521-2324  
Email: [stephanie.greeman@roche.com](mailto:stephanie.greeman@roche.com)

Date Prepared: September 23, 2010

---

**Device Name** Proprietary names: Elecsys<sup>®</sup> SHBG Immunoassay System

Common name: SHBG test

Classification name: Radioimmunoassay, Testosterones and Dihydrotestosterone

---

*Continued on next page*

## 510(k) Summary, continued

---

**Device  
Description**

A device for the measurement of human SHBG in serum or plasma.

---

**Substantial  
Equivalence**

The Elecsys SHBG Test System is substantially equivalent to other devices legally marketed in the United States. We claim equivalency to the currently marketed Elecsys SHBG Test System (K031717).

---

*Continued on next page*

## 510(k) Summary, continued

<b>Premarket Notification, Traditional 510(k) for Elecsys SHBG Assay</b>		
<b>Immunoassay Comparison</b>		
<b>Feature</b>	<b>Predicate Device: Elecsys SHBG Assay (K031717)</b>	<b>Elecsys SHBG Assay</b>
<b>General Assay Features</b>		
<b>Intended Use/ Indications for Use</b>	Immunoassay for the <i>in vitro</i> quantitative determination of sex hormone-binding globulin in human serum and plasma.  The electrochemiluminescence immunoassay "ECLIA" is intended for use on Elecsys and <b>cobas e</b> immunoassay analyzers.	Same
<b>Assay Protocol</b>	Sandwich assay	Same
<b>Detection Protocol</b>	Electrochemiluminescent Immunoassay	Same
<b>Application</b>	18 Minute	Same
<b>Instrument Platform</b>	Roche Elecsys 1010, 2010/ <b>cobas e</b> 411 and MODULAR ANALYTICS E170/ <b>cobas e</b> 601	Same with the exception of the removal of the Elecsys 1010 analyzer
<b>Sample Volume</b>	10 $\mu$ L	Same
<b>Sample Type</b>	Human serum and plasma treated with lithium heparin.	Same
<b>Traceability</b>	The Elecsys SHBG assay has been standardized against the 1 <sup>st</sup> International Standard for SHBG, NIBSC code 95/560.	Same
<b>Calibrator</b>	Elecsys SHBG CalSet	Same

*Continued on next page*

## 510(k) Summary, continued

Immunoassay Comparison		
Feature	Predicate Device: Elecsys SHGB Assay (K031717)	Elecsys SHBG Assay
<b>General Assay Features</b>		
<b>Calibration Interval</b>	<p>Calibration must be performed once per reagent lot using fresh reagent (i.e. not more than 24 hours since the reagent kit was registered on the analyzer). Renewed calibration is recommended as follows:</p> <p>MODULAR ANALYTICS E170, Elecsys 2010 and <b>cobas e</b> analyzers:</p> <ul style="list-style-type: none"> <li>• After 1 month (28 days) when using the same reagent lot</li> <li>• After 7 days (when using the same reagent kit on the analyzer)</li> </ul> <p>Elecsys 1010 analyzer:</p> <ul style="list-style-type: none"> <li>• With every reagent kit</li> <li>• After 7 days (ambient temperature 20-25 °C)</li> <li>• After 3 days (ambient temperature 25-32 °C)</li> </ul>	Same with the exception of the removal of the Elecsys 1010 analyzer.
<b>Controls</b>	Elecsys PreciControl Universal 1 and 2	Same

*Continued on next page*

## 510(k) Summary, continued

Immunoassay Comparison		
Feature	Predicate Device: Elecsys SHBG Assay (K031717)	Elecsys SHBG Assay
<b>General assay features</b>		
<b>Reagent Stability/Storage</b>	<p>Store at 2-8 °C. Store the Elecsys SHBG reagent kit upright in order to ensure complete availability of the microparticles during automatic mixing prior to use.</p> <p>Stability:            Unopened at 2-8 °C—up to the stated expiration date            After opening at 2-8 °C—12 weeks            On MODULAR ANALYTICS E170 and <b>cobas e 601</b>—7 weeks            On Elecsys2010 and <b>cobas e 411</b>—7 weeks            On Elecsys 1010- 4 weeks (stored alternately in the refrigerator and on the analyzer – ambient temperature 20-25 °C; up to 20 hours opened in total)</p>	Same with the exception of the removal of the Elecsys 1010 analyzer.

*Continued on next page*

## 510(k) Summary, continued

Immunoassay Comparison						
Feature	Predicate Device: Elecsys SHBG Assay (K031717)	Elecsys SHBG Assay				
Labeled Performance Characteristics						
<b>Measuring Range</b>	0.350-200 nmol/L	0.800 -200 nmol/L				
<b>Expected Values</b>	Males: 10-80 nmol/L Females, non-pregnant: 20-130 nmol/L	Reference ranges for the SHBG (nmol/L)				
			N	Median	5 <sup>th</sup> Perc	95 <sup>th</sup> perc
		Males 20-49 years	136	33.5	16.5	55.9
		Males ≥ 50 years	78	40.8	19.3	76.4
		Females 21-49 years	89	64.3	24.6	122
Females ≥ 50 years	71	57.4	17.3	125		

*Continued on next page*

# 510(k) Summary, continued

Immunoassay Comparison					
Feature	Predicate Device: Elecsys SHBG Assay (K031717) (18 Minute)	Elecsys SHBG Assay			
Labeled Performance Characteristics					
Expected Values cont.	Reference ranges for Testosterone II (nmol/L)				
		N	Median	5 <sup>th</sup> Perc	95 <sup>th</sup> perc
	Males 20-49 years	136	18.6	8.64	29.0
	Males ≥ 50 years	78	16.5	6.68	25.7
	Females 21-49 years	89	0.941	0.290	1.67
	Females ≥ 50 years	71	0.563	0.101	1.42
	Calculation for obtaining FTI (or FAI): %FTI = (Testosterone (nmol/L) ÷ SHBG (nmol/L)) x 100				
	Reference ranges for Free Testosterone Index (FTI) / Free Androgen Index (FAI) (nmol/L)				
		N	Median	5 <sup>th</sup> Perc	95 <sup>th</sup> perc
	Males 20-49 years	136	57.2	35.0	92.6
Males ≥ 50 years	78	38.2	24.3	72.1	
Females 21-49 years	89	1.53	0.297	5.62	
Females ≥ 50 years	71	1.15	0.187	3.63	

Continued on next page

<b>Immunoassay Comparison</b>		
<b>Feature</b>	<b>Predicate Device: Elecsys SHBG Assay (K031717) (18 Minute)</b>	<b>Elecsys SHBG Assay</b>
<b>Labeled Performance Characteristics</b>		
<b>Precision</b>	E1010/ 2010/ 411 <ul style="list-style-type: none"> <li>• Within Run 2.1 – 2.7 % CV from 14.1 -204 nmol/L</li> <li>• Total 2.6 -5.6% CV from 14.1 -204 nmol/L</li> </ul> E170/ e601 <ul style="list-style-type: none"> <li>• Within run 1.1 – 1.7% CV from 14.9 – 219 nmol/L</li> <li>• Total 1.8 – 4.0% CV from 14.9 -219 nmol/L</li> </ul>	Same with the exception of the removal of the Elecsys 1010 analyzer.
<b>Analytical Sensitivity</b>	Limit Detection Level (LDL): 0.35 nmol/L	Limit of Blank (LoB): 0.500 nmol/L  Limit of Detection (LoD): 0.800 nmol/L  Limit of Quantitation (LoQ): 2.00 nmol/L
<b>Hook Effect</b>	There is no high-dose hook effect at SHBG concentrations up to 1000 nmol/L.	Same

---

*Continued on next page*

## 510(k) Summary, continued

Immunoassay Comparison		
Feature	Predicate Device: Elecsys SHBG Assay (K031717) (18 Minute)	Elecsys SHBG Assay
Labeled Performance Characteristics		
<b>Limitations</b>	<p>The assay is unaffected by:</p> <ul style="list-style-type: none"> <li>• Hemoglobin &lt; 2.9 mg/dL</li> <li>• Bilirubin &lt; 60 g/dL</li> <li>• Intralipid &lt; 2700 mg/dL</li> <li>• Biotin &lt; 60 ng/mL</li> <li>• Rheumatoid factors up to 1160 IU/mL</li> <li>• In vitro tests were performed on 16 commonly used pharmaceuticals. No interference with the assay was found.</li> <li>• As with all tests containing monoclonal mouse antibodies, erroneous findings may be obtained from samples taken from patients who have been treated with monoclonal mouse antibodies or have received them for diagnostic purposes.</li> <li>• In rare cases, interference due to extremely high titers of antibodies to streptavidin and ruthenium can occur. These effects are minimized by suitable test design.</li> <li>• For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.</li> </ul>	Same

*Continued on next page*

## 510(k) Summary, continued

Immunoassay Comparison			
Feature	Predicate Device: Elecsys SHBG Assay (K031717) (18 Minute)	Elecsys SHBG Assay	
<b>Labeled Performance Characteristics</b>			
<b>Method Comparison</b>	A comparison of the Elecsys SHBG assay (y) with a commercially available SHBG assay (x) using clinical samples gave the following correlation:		
	n = 109	<b>Passing/Bablok</b>	<b>Linear Regression</b>
	Min = 11.2 nmol/L		
	Max = 155 nmol/L		
	Slope	1.17	1.15
	Intercept	-3.26	-1.82
	Tau / r	0.909	0.981

---

**Confidentiality** Roche Diagnostics Corporation requests that the FDA not disclose the nature or existence of the premarket notification until the substantial equivalence decision has been reached.

---

**Closing** We trust that the information provided in this Premarket Notification will support a determination of substantial equivalence for the Elecsys SHBG Immunoassay.

If you should have questions or require further information, please do not hesitate to contact me.

- Phone: (317) 521-3208
- FAX: (317) 521-2324
- [kelly.french@roche.com](mailto:kelly.french@roche.com)

Sincerely,

Kelly French, RN, BSN, RAC  
 Regulatory Affairs Consultant  
 US Regulatory Submissions  
 Roche Diagnostics Corporation



Food and Drug Administration  
10903 New Hampshire Avenue  
Silver Spring, MD 20993

Roche Diagnostics  
c/o Ms. Kelly French  
Regulatory Affairs Consultant  
9115 Hague Road  
P.O. Box 50416  
Indianapolis, IN 46250-0416

JAN 20 2011

Re: k102814  
Trade Name: Elecsys SHBG  
Regulation Number: 21 CFR §862.1680  
Regulation Name: Radioimmunoassay, testosterone and dihydrotestosterone  
Regulatory Class: Class I, reserved  
Product Codes: CDZ  
Dated: December 15, 2010  
Received: December 16, 2010

Dear Ms. French:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

If you desire specific advice for your device on our labeling regulation (21 CFR Parts 801 and 809), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read 'CHC', followed by a long horizontal line extending to the right.

Courtney Harper, Ph.D.  
Director  
Division of Chemistry and Toxicology  
Office of *In Vitro* Diagnostic Device  
Evaluation and Safety  
Center for Devices and Radiological Health

Enclosure

## Indication for Use

510(k) Number (if known): K102814

Device Name: Elecsys SHBG

Indication For Use:

Immunoassay for the in vitro quantitative determination of sex hormone-binding globulin in human serum and plasma. The Elecsys SHBG Immunoassay is intended for use as an aid in the diagnosis of androgen disorders. The electrochemiluminescence immunoassay "ECLIA" is intended for use on Elecsys and **cobas e** immunoassay analyzers.

Prescription Use  X   
(21 CFR Part 801 Subpart D)

And/Or

Over the Counter Use        
(21 CFR Part 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE; CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Device Evaluation and Safety (OIVD)



Division Sign-Off  
Office of In Vitro Diagnostic Device  
Evaluation and Safety

510(k)  K102814