



K 103 206

MAR - 4 2011

RAUMEDIC AG · Postfach 5 01 · D - 95205 Münchberg

**Forschung & Entwicklung**

Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center – W066-0609  
10903 New Hampshire Avenue  
Silver Spring, MD 20993-0002

Telefon: (0 92 52) 359-0  
Telefax: (0 92 52) 359-10 00  
info@RAUMEDIC.com  
www.RAUMEDIC.de

Direktkontakt:  
Telefon: (0 92 52) 359-2782

Reiner.Thiem@RAUMEDIC.com

Ihre Nachricht vom  
Ihre Zeichen

Sachbearbeiter H Thiem  
Unsere Zeichen thi

Tag  
22.02.2011

**510(k) Summary**

**This summary of 510(k) safety and effectiveness information is submitted in accordance with the requirements of SMDA 1990 and 21 CFR §807.92.**

**1. Name, address, phone and fax number of the applicant**

RAUMEDIC AG  
Hermann-Staudinger-Straße 2  
95233 Helmbrechts  
D - Germany  
Tel.: 0049/9252/359-0  
Fax: 0049/9252/359-1000

**2. Contact person**

Mr. Reiner Thiem  
Head of Regulatory Affairs  
Hermann-Staudinger-Straße 2  
95233 Helmbrechts  
D – Germany  
Tel.: 0049/9252/359-2782

**3. Date of preparation of the summary**

August, the 26<sup>th</sup> 2010



#### 4. Name of the device

The **RAUMEDIC® -ICP-Monitoring-System** is composed of the following elements:

- Neurovent® - P
- Neurovent® - P - C
- RAUMEDIC® - Bolt Kit (CH 5)
- RAUMEDIC® - Drill Kit 4,5 mm (CH 5)
- ICP-TEMP-Cable
- NPS2
- NPS3

|                             |   |
|-----------------------------|---|
| Device Classification Name: | Device, Monitoring, Intracranial Pressure |
| Classification Panel:       | Neurology                                 |
| CFR Section:                | 21 CFR §882.1620                          |
| Device Class:               | Class II                                  |
| Product Code:               | GWM                                       |

#### 5. Device Description

The **RAUMEDIC® -ICP-Monitoring-System** determines safely, quickly and accurately the level and change in intracranial pressure (ICP) by using semi-conductor pressure sensors.

The Neurovent® - P and Neurovent® - P - C are indicated for use in parenchymal pressure monitoring. Both types of catheters are implanted in parenchyma via a RAUMEDIC® - Bolt Kit CH5. In addition to the catheter a zero point module NPS2 x is needed. "x" depends on the type of patient monitor available in the hospital - there are 20 different references. To the equipment also belongs an ICP-Temp-Cable and a NPS3 pressure display unit. The difference between the Neurovent® - P and the Neurovent® - P - C is that the housing material in the C-version is ceramic instead of titanium.

The **RAUMEDIC® -ICP-Monitoring-System** is composed of the following elements:

- Neurovent® - P
- Neurovent® - P - C
- RAUMEDIC® - Bolt Kit (CH 5)
- RAUMEDIC® - Drill Kit 4,5 mm (CH 5)
- ICP-TEMP-Cable
- NPS2
- NPS3

#### 6. Device Intended Use

The **RAUMEDIC® -ICP-Monitoring-System** is indicated for use in parenchymal pressure monitoring and can be used for the measurement of the intra-cranial pressure (ICP) as well as of the cerebral perfusion pressure (central arterial blood pressure minus ICP) which is the essential pre-requisite for an effective treatment of patients suspected of suffering from intracranial pressure increases (such as cranio-cerebral traumas, GCS  $\leq$  8; malignant medial cardiac infarctions; hepatic encephalopathy; SAB Hunt / Hess IV + V; cerebral edema; hydrocephalus) or of patients whose clinical picture may be linked to an increase of the ICP.



Use of the parenchymal intracranial pressure monitoring kit with bolt is contra-indicated in children under one year old.

The RAUMEDIC® precision pressure catheters are MR Unsafe.

## 7. Substantial Equivalence Summary

The RAUMEDIC® -ICP-Monitoring-System is substantially equivalent to those of a legally marketed predicate device, the Pressio® ICP MONITORING SYSTEM, which was cleared to market under 510 (k) K062584 on July 5.  
See device comparison table attached.

Based on performance testing and the available information concerning the referenced comparison device, the **RAUMEDIC® -ICP-Monitoring-System** is similar in that:

- The devices have the same intended use and indication for use.
- The devices are made of the same materials or substantially similar materials.
- The devices have similar form, function, procedures and features.
- Performance characteristics are suitable for designated indications for use

Based on this, the anticipated clinical performance of the **RAUMEDIC® -ICP-Monitoring-System** is equivalent to the referenced systems.

## 8. Device Testing

Biocompatibility studies were conducted per ISO 10993 standard and have demonstrated that the materials used to manufacture the RAUMEDIC® -ICP-Monitoring-System are safe for its intended use.

In addition, the two implantable catheters were subjected to extensive performance testing. Results of the testing showed that the catheter designs are safe for their intended uses.

The NPS3 underwent numerous safety tests, including testing to IEC 60601-1.

Finally, the manufacturing process of the RAUMEDIC® - ICP- Monitoring -System complies with the United States Food and Drug Administration and European Standards for the manufacturing of medical devices.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - WO66-G609  
Silver Spring, MD 20993-0002

Raumedic AG  
c/o Mr. Norbert Stuibler  
TUV SUD Product Service GmbH  
1775 Old Highway 8, NW  
Suite 104  
New Brighton, MN 55112-1891

MAR - 4 2011

Re: K103206

Trade/Device Name: RAUMEDIC® - ICP Monitoring System  
Regulation Number: 21 CFR 882.1620  
Regulation Name: Intracranial Pressure Monitoring Device  
Regulatory Class: Class II  
Product Code: GWM  
Dated: February 25, 2011  
Received: February 28, 2011

Dear Mr. Stuibler:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Malvina B. Eydelman, M.D.  
Director  
Division of Ophthalmic, Neurological,  
and Ear, Nose and Throat Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

K103206

### Indications for Use

510(k) Number (if known): K \_\_\_\_\_

Device Name: Device, Monitoring, Intracranial Pressure

Indications for Use:


The RAUMEDIC® -ICP-Monitoring-System is indicated for use in parenchymal pressure monitoring.

Use of the parenchymal intracranial pressure monitoring kit with bolt is contra-indicated in children under one year old.

The RAUMEDIC® precision pressure catheters are MR Unsafe.

Prescription Use  AND/OR Over-The-Counter Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)



(Division Sign-Off)  
Division of Ophthalmic, Neurological and Ear,  
Nose and Throat Devices

Prescription Use  \_\_\_\_\_  
(Per 21 CFR 801.109)

510(k) Number K103206

---