

**Attachment II 510(k) Summary**

As Required by CFR 870.92

<b>Preparing Date</b>	November 23, 2010	FEB - 4 2011
<b>Sponsor</b>	Shenzhen Jumper Medical Equipment Co., Ltd 5th Floor, Building No.34, Baoyuan Industrial Zone, Xixiang Street, Baoan District, Shenzhen, Guangdong 518067, P.R. China	
	Boley He      Quality Manager Tel: +86-755-2669 6279 Fax: +86-755-2685 2678 Email: <a href="mailto:rd@jumper-medical.com">rd@jumper-medical.com</a>	
<b>Submission Correspondent</b>	Ms. Diana Hong / Mr. Lee Fu Shanghai Mid-Link Business Consulting Co., Ltd P.O. BOX 237-023 Shanghai, 200237, China Tel: +86-21-64264467 Fax: 240-238-7587 Email: <a href="mailto:diana.hong@mid-link.net">diana.hong@mid-link.net</a>	
<b>Proposed Product</b>		
Trade Name	Fetal Doppler	
Model	JPD-100S	
Product Code:	KNG	
Regulation Number:	21 CFR 884.2660	
Device Class:	Class II	
<b>Submission Purpose:</b>	New Device	
<b>Submission Type</b>	Traditional 510(k)	
<b>Predicate Device:</b>		
Device Name	Contec Pocket Fetal Doppler Baby Sound A	
K Number	K082480	
Manufacturer	Contec Medical Systems Co., LTD	
<b>Intended Use</b>	The Fetal Doppler JPD-100S is a hand-held, battery powered audio	

<b>Device Description</b>	<p>Doppler device used for detecting fetal heart beats.</p> <p>JPD-100S Fetal Doppler is a prescription from licensed physician in hospitals, clinics and private offices. It is a hand-held, battery powered audio Doppler device integrated with 3MHz probe, used for detecting fetal heart beats. And the device is normally applied to pregnancy for more than twelve weeks.</p>
<b>Test Conclusion</b>	<p>Laboratory testing was conducted to validate and verify that the proposed devices met all design specifications, including electrical safety, EMC, specification per following standards.</p> <p>IEC60601-1(1988; Amendment 1, 1991-11, Amendment 2, 1995): Medical electrical equipment –Part 1: General requirements for basic safety and essential performance</p> <p>IEC60601-1-2 (2001): Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility -Requirements and tests</p> <p>IEC 61266 (1994 ): Ultrasonic - Hand-held probe Doppler fetal heartbeat detectors -Performance requirements and methods of measurement and reporting</p> <p>ISO 10993-5 (2009): Biological evaluation of medical devices -- Part 5: Tests for In Vitro cytotoxicity</p> <p>ISO 10993-10 (2002; Amd. 1:2006): Biological evaluation of medical devices - Part 10: Tests for irritation and delayed-type hypersensitivity</p> <p>Results of these tests demonstrate compliance to the requirements of all consensus standards.</p>
<b>Comparison Summary</b>	<p>The proposed device shares same classification, intended use, similar technical specifications, same safety specifications and skin-contact material with the predicate device. The main difference includes battery, operation and storage conditions, transducer frequency and /sata, but the difference is considered to have no effect on effectiveness and safe. The SE analysis of the difference is provided as following.</p> <p>The batteries of proposed and predicate device are different. But both of them comply with the requirements of electrical safety and EMC of IEC 60601-1 and IEC 60601-1-2. Therefore, this difference is</p>

considered has no effect on effectiveness and safety.

The operation and storage conditions of the proposed device have slight differences with those of the predicate device. But these conditions are considered to be able to cover the normal working and storage environment. In addition, these conditions are clearly on the proposed labeling, to remind the user handle and operate with these conditions. Therefore, this difference is considered to have no effect on effectiveness and safe.

The transducer frequency of proposed device and predicate device is different. But the proposed transducer was tested for its acoustic output, and the result demonstrates that it complies with the standard criteria. In addition, 3MHz transducer of fetal doppler is very common in clinical use, such as, 3MHz normal transducer of Contec Pocket Fetal Doppler SONOLINE A and SONOLINE B (k082480). Therefore, this difference is considered to have no effect on effectiveness and safe.

The  $I_{sata}$  of proposed device and predicate device is different. But the proposed transducer was tested for its acoustic output, and the result demonstrates that it complies with the standard criteria. In addition, transducer of fetal doppler with  $I_{sata} < 5 \text{ mW/cm}^2$  is very common in clinical use, such as, CMS-150-T2 and CMS-150-T5 transducer of Contec Pocket Fetal Doppler SONOLINE A and SONOLINE B (k082480). Therefore, this difference is considered to have no effect on effectiveness and safe.

**SE Determination**

The proposed device, Fetal Doppler JPD-100S, is substantially equivalent (SE) to the predicate device Contec Pocket Fetal Doppler Baby Sound A (K082480).



Food and Drug Administration  
10903 New Hampshire Avenue  
Silver Spring, MD 20993

Shenzhen Jumper Medical Equipment Co., Ltd  
% Mr. Marc M. Mouser  
Engineering Leader/FDA Office Coordinator  
Underwriters Laboratories, Inc.  
2600 N.W. Lake Road  
CAMAS WA 98607-8542

FEB - 4 2011

Re: K110124

Trade/Device Name: Fetal Doppler JPD-100s  
Regulation Number: 21 CFR 884.2660  
Regulation Name: Fetal ultrasound monitor and accessories  
Regulatory Class: II  
Product Code: KNG  
Dated: December 17, 2010  
Received: January 18, 2011

Dear Mr. Mouser:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the Fetal Doppler JPD-100s, as described in your premarket notification:

Transducer Model Number

3.0MHz CW

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucml15809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

If you have any questions regarding the content of this letter, please contact Lauren Hefner at (301) 796-6881.

Sincerely Yours,



Mary Pastel, ScD.  
Director  
Division of Radiological Devices  
Office of In Vitro Diagnostic Device  
Evaluation and Safety  
Center for Devices and Radiological Health

Enclosure(s)

## Section II Indications for Use Statement

510(k) Number: \_\_\_\_\_

Device Name: \_\_\_\_\_

Fetal Doppler JPD-100S

### Indications for Use:

The Fetal Doppler JPD-100S is a hand-held, battery powered audio Doppler device used for detecting fetal heart beats.

Prescription Use  \_\_\_\_\_  
(Part 21 CFR 801 Subpart D)

AND/OR Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THE LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)



(Division Sign-Off)

Division of Radiological Devices

Office of *In Vitro* Diagnostic Device Evaluation and Safety

510(k) Number K110124

Page 1 of 2

**Diagnostic Ultrasound Indications for Use Form**

**Device Name:** Fetal Doppler JPD-100S

**Transducer:** 3.0MHz CW inherent of the main unit


**Intended Use:** Detect fetal heart beats as follows:

Clinical Application		Mode of Operation						
General (Track 1 only)	Specific (Track 1&3)	B	M	PWD	CWD	Color Doppler	Combined (Specify)	Other* (Specify)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal				N			
	Abdominal							
	Intra-operative (Specify)							
	Intra-operative (Neuro)							
	Laparoscopic							
	Pediatric							
	Small Organ (Specify)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skeletal (Conventional)							
	Musculo-skeletal (Superficial)							
Intravascular								
Other (Specify)								
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Intravascular (Cardiac)							
	Trans-esoph. (Cardiac)							
	Intra-cardiac							
Other (Specify)								
Peripheral Vessel	Peripheral vessel							
	Other (Specify)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E Comments:

**(PLEASE DO NOT WRITE BELOW THE LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)**

Concurrence of CDRII, Office of Device Evaluation (ODE)

  
(Division Sign-Off)  
Division of Radiological Devices  
Office of In Vitro Diagnostic Device Evaluation and Safety  
510K K110621