

**510(k) Summary**  
per 21 CFR §807.92 (c)

MAR - 3 2011

K110295

|                                     |   |         |                 |
|-------------------------------------|---|---------|-----------------|
| <b>Submitter's Name and Address</b> | Boston Scientific Corporation<br>One Scimed Place<br>Maple Grove, MN 55311  |         |                 |
| <b>Contact Name and Information</b> | Shannon Pettit<br>Senior Regulatory Affairs Specialist<br>Tel: 763-494-2833<br>Fax: 763-494-2222<br>E-mail: Shannon.Pettit@bsci.com   |         |                 |
| <b>Date Prepared</b>                | January 31, 2011  |         |                 |
| <b>Trade Name</b>                   | Interlock™-35 Fibered IDC™ Occlusion System   |         |                 |
| <b>Common Name</b>                  | Vascular embolization device  |         |                 |
| <b>Classification</b>               | Class II  |         |                 |
| <b>Product Code</b>                 | KRD, Vascular embolization devices<br>(21 CFR 870.3300)   |         |                 |
| <b>Predicate Devices</b>            | Boston Scientific   | K060078 | SE: 31 Jan 2006 |
|                                     | Interlock Fibered IDC Occlusion System  |         |                 |
| <b>Predicate Devices</b>            | Fibered Platinum  | K955293 | SE: 6 Feb 1996  |
|                                     | .035" Type Occlusion Coils  |         |                 |
| <b>Reason for Submission</b>        | To gain clearance for a modified fibered interlocking detachable coil embolization system based on the currently marketed Fibered IDC Occlusion System cleared under K060078. |         |                 |

---

**Device  
Description**

The Interlock-35 Fibered IDC Occlusion System includes a coil manufactured from a platinum-tungsten alloy that is mechanically attached to a coil delivery wire. This assembly is contained within an introducer sheath. The platinum coil contains synthetic fibers for greater thrombogenicity. The Interlock-35 Fibered IDC Occlusion System is designed to be delivered under fluoroscopy through a 5F (0.035 in [0.89 mm] or 0.038 in [0.97 mm] inner lumen) Diagnostic Catheter (e.g. Imager™ II Diagnostic Catheter). The interlocking delivery wire design allows the coil to be advanced and retracted before final placement in the vessel, thus aiding in more controlled delivery including the ability to withdraw the coil prior to deployment.

---

**Indications for  
Use**

The Interlock-35 Fibered IDC Occlusion System is indicated to obstruct or reduce rate of blood flow in the peripheral vasculature. This device is not intended for neurovascular use.

---

**Comparison of  
Technological  
Characteristics**

The proposed Interlock-35 Fibered IDC Occlusion System is a modification to existing designs/shapes currently marketed by Boston Scientific. It is important to note that each of the Interlock-35 Fibered IDC Occlusion System device characteristics is currently available on BSC legally marketed devices and no changes in fundamental scientific technology have been made. Additionally, the performance characteristics of the Interlock-35 Fibered IDC Occlusion System are substantially equivalent to the currently marketed devices.

---

**Non-Clinical  
Performance  
Data**

---

Determination of substantial equivalence is based on an assessment of non-clinical performance data. Non-clinical performance data submitted in support of substantial equivalence is based on the Failure Modes/Effects Analysis (FMEA) risk analysis method completed for the Interlock-35 Fibered IDC Occlusion System to demonstrate that the proposed devices are suitable for their intended use.

All testing performed and data demonstrate passing results according to executed verification protocols. Therefore, results of non-clinical performance data submitted supports substantial equivalence to predicate device.

The following performance tests were performed:

- Coil Geometry
- Stretch Resistance
- Deployment
- Atraumatic Embolic Coil Tip
- Fiber Retention
- Interlocking Arm to Embolic Coil Tensile Strength
- Embolic Coil Corrosion Resistance
- Radial Force
- Atraumatic Delivery Wire Tip
- Interlocking Arm to Delivery Wire Tensile Strength
- Corrosion Resistance
- Delivery Wire Coating Fracture
- Delivery Wire Coating Flexure
- Delivery Wire OD
- Delivery Wire Length
- Introducer Sheath Interface
- Introducer Sheath Deliverability
- Delivery Catheter Compatibility
- Rotating Hemostatic Valve (RHV) Compatibility

**Clinical  
Performance  
Data**

---

Not applicable; determination of substantial equivalence is not based on an assessment of clinical performance data.

---

**Conclusion**

---

The modifications do not affect the intended use or alter the fundamental scientific technology of the predicate Boston Scientific Interlock Fibered IDC Occlusion System (K060078).

Based on the Indications for Use, unaltered technological characteristics, and submitted non-clinical performance data supporting this modification, the Boston Scientific Interlock-35 Fibered IDC Occlusion System is shown to be appropriate for its intended use and demonstrates that the device is as safe, as effective, and performs as well as the predicate device.

---



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - WO66-G609  
Silver Spring, MD 20993-0002

MAR - 3 2011

Boston Scientific Corporation  
c/o Ms. Shannon Pettit  
Senior Regulatory Affairs Specialist  
One Scimed Place  
Maple Grove, MN 55311

Re: K110295  
Trade/Device Name: Interlock™-35 Fibered IDC™ Occlusion System  
Regulation Number: 21 CFR 870.3300  
Regulation Name: Vascular Embolization Device  
Regulatory Class: Class II  
Product Code: KR D  
Dated: January 31, 2011  
Received: February 1, 2011

Dear Ms. Pettit:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be

Page 2 - Ms. Shannon Pettit


found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



 Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K110295

Device Name: Interlock-35™ Fibered IDC™ Occlusion System

### Indications for Use:

The Interlock-35 Fibered IDC Occlusion System is indicated to obstruct or reduce rate of blood flow in the peripheral vasculature. This device is not intended for neurovascular use.

Prescription Use X AND/OR Over-The-Counter Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of 1

Diana R. Vachon  
(Division Sign-Off)  
Division of Cardiovascular Devices

510(k) Number K110295