



K110626

APR 20 2011

## 510(k) SUMMARY

### General Information:

Date of Summary Preparation: March 3, 2011

Name and Address of Manufacturer: Pathway Medical Technologies, Inc.  
10801 120<sup>th</sup> Ave NE  
Kirkland, Washington 98033

Contact Person: Brit Baird  
Regulatory Affairs Manager  
Phone: 425-636-4137  
Fax: 425-636-4001

Device Trade Names: JETSTREAM Navitus™ System  
JETSTREAM G3® SF System  
JETSTREAM G3® SE System  
JETSTREAM G3® L System

Common Name: Peripheral Atherectomy Catheter

Regulation Number: 21 CFR 870.4875

Regulation Name: Intraluminal Artery Stripper

Regulatory Class: Class II

Classification Panel: Cardiovascular

Product Code: MCW

000066

**Performance Standards:** Performance Standards do not currently exist for these devices. None are established under Section 514.

**Device Description:** The Jetstream Systems are atherectomy catheter systems designed with either a fixed (Jetstream G3 SF System) or an expandable (Jetstream Navitus, G3 SE and G3 L Systems) cutting tip intended for use in debulking and treating vascular disease in the peripheral vasculature. Separate lumens within the Catheter allow for continuous aspiration and infusion during device use. Excised tissue, thrombus, and fluid are aspirated from the peripheral treatment site through a port in the Catheter tip to an external collection bag located on the Console. The distal portion of the Catheter also possesses infusion ports that provide continuous infusion of sterile saline during the atherectomy procedure.

The Jetstream Systems consist of two primary components: a Catheter with Control Pod and a Console, which are packaged separately. Each of these system components is described generally as follows:

- **Jetstream Catheter with Control Pod:** A sterile, single-use unit consisting of an electrically-driven Catheter with attached Control Pod. As with the predicate devices, the Jetstream Catheters continue to utilize a differentially cutting tip and include both aspiration and infusion capabilities, and the Control Pod provides a user interface with keypad controls. The unit, its electrical connectors, tubing, and aspirant collection bag are packaged in a double pouched tray.
- **PV Console:** A reusable compact PV Console, with two (2) peristaltic pumps for aspiration and infusion, power supply, system controller, keypad interface, and LED indicators for device operational status. The PV Console mounts on a standard I.V. stand and remains outside the sterile field during the procedure.

The primary modifications of this 510(k) simplify the user interface by eliminating the two buttons that increase and decrease speed on the Control Pod Membrane. These modifications apply to the entire family of Jetstream Systems. In addition, non-significant changes made to the Jetstream G3 System (which is being given a new trade name - "Jetstream Navitus System" - with the modifications included in this 510(k)) and PV Console are also included.

**Indications for Use:** The JETSTREAM System is intended for use in atherectomy of the peripheral vasculature and to break apart and remove thrombus from upper and lower extremity peripheral arteries. It is not intended for use in coronary, carotid, iliac or renal vasculature.

000067

**Substantially Equivalent Devices:** Pathway Medical cites the following devices as substantially equivalent predicate devices:

Primary Predicate Device	Pathway Medical Predicate 510(k)
JETSTREAM G3® System	K101221
JETSTREAM G3® SF System	K101334
JETSTREAM G3® SE System	K101334
JETSTREAM G3® L System	K100462
JETSTREAM G3® L System	K093918
JETSTREAM G2™ NXT System	K091509
JETSTREAM Pathway PV™ Atherectomy System	K082186
Pathway PV™ Atherectomy System	K081328

**Testing Summary:** To demonstrate substantial equivalence of the modified Jetstream Systems to the predicate devices, the technological and performance characteristics were evaluated using *in vitro* testing for the primary (and supporting) modifications, as outlined below:

- System Reliability/Life Test
- Main Cable Joint Strength
- Control Pod Electrical Testing
- Control Pod Logic Verification
- Electrical Safety
- Intended Use

In addition, the following *in vitro* testing was provided to support the previously mentioned non-significant changes to the predicate Jetstream G3 System:

- System Reliability/Life Test
- Accessory Compatibility
- Infusion & Aspiration Flow Rates
- Aspiration Efficiency & Crossing Time
- Rotational Speed
- Material Liberation (Teflon & Polyimide)
- Aggressive Thrombus Testing

The results from these tests:

- demonstrate that the technological and performance characteristics of the modified Jetstream Systems are comparable to the predicate devices,
- support the safety and effectiveness of the modifications that are the subject of this 510(k), and
- ensure the modified device can perform in a manner equivalent to the listed predicate devices with the identical intended use.

000068

**Conclusion (Statement of Equivalence):** The data and information presented within this submission (including *in vitro* testing) and the similarities between the modified and predicate devices support a determination of substantial equivalence, and therefore market clearance of the modified Jetstream Systems through this 510(k) Premarket Notification.

000069



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room -WO66-G609  
Silver Spring, MD 20993-0002

Pathway Medical Technologies, Inc.  
c/o Mr. Brian Cleary  
Vice President of Clinical and Regulatory Affairs  
10801 120<sup>th</sup> Ave NE  
Kirkland, WA 98033

APR 20 2011

Re: K110626

Trade/Device Name: Jetstream Navitus, G3 SF, G3 SE and G3 L Systems

Regulation Number: 21 CFR 870.4875

Regulation Name: Intraluminal Artery Stripper

Regulatory Class: Class II (two)

Product Code: MCW

Dated: April 7, 2011

Received: April 8, 2011

Dear Mr. Cleary:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucml15809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

INDICATIONS FOR USE

510(k) Number (if known): K110626


Device Name: JETSTREAM Navitus™ System  
JETSTREAM G3® SF System  
JETSTREAM G3® SE System  
JETSTREAM G3® L System

**Indications for Use:** The JETSTREAM System is intended for use in atherectomy of the peripheral vasculature and to break apart and remove thrombus from upper and lower extremity peripheral arteries. It is not intended for use in coronary, carotid, iliac or renal vasculature.

Prescription Use  X  AND/OR Over-The-Counter Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
\_\_\_\_\_  
(Division Sign-Off)  
Division of Cardiovascular Devices  
510(k) Number K110626