

K110723

## 510(k) SUMMARY

JUL 28 2011

**Submitted by:** Masimo Corporation  
40 Parker, Irvine, CA 92618  
Phone: 949-297-7000  
Fax: 949-297-7592

**Company Contact:** Anil Bhalani, Director of Regulatory Affairs

**Date Summary Prepared:** May 11, 2011

**Trade Name** M-LNCS and Reprocessed M-LNCS Oximetry Sensors

**Regulation Number:** 21 CFR 870.2700

**Regulation Class:** Class II

**Regulation Name:** Oximeter

**Common Name:** Oximeter Sensor

**Product Code:** NLF, DQA, DSA

**Substantially Equivalent Devices:** Masimo LNCS/M-LNCS Oximetry Sensors, 510(k) No. K101896  
LNCS Oximetry Sensors-Sterile, 510(k) No. K083622  
Reprocessed LNCS Oximetry Sensors, 510(k) No. K083719

### Device Description

The M-LNCS and Reprocessed M-LNCS Oximetry Sensors are fully compatible disposable sensors for use with instruments which include or compatible with the following technologies:

- Masimo SET technology
- Masimo Rainbow SET technology
- Nellcor technology

The M-LNCS series has been validated with Masimo SET Oximetry Technology and on Nellcor's N-200 Pulse Oximeter. The saturation accuracy of the Neonate and Preterm sensors were validated on adult volunteers and 1% was added to account for properties of fetal hemoglobin.

### Intended Use/ Indications for Use

The M-LNCS Sensors are indicated for the continuous non-invasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO<sub>2</sub>) and pulse rate (measured by an SpO<sub>2</sub> sensor) for use with adult, pediatric, infant, and neonatal patients during both no motion and motion conditions, and for patients who are well or poorly perfused in hospitals, hospital-type facilities, mobile, and home environments.

The Reprocessed M-LNCS Sensors are indicated for the continuous non-invasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO<sub>2</sub>) and pulse rate (measured by an SpO<sub>2</sub> sensor) for use with adult, pediatric, infant, and neonatal patients during both no motion and motion conditions, and for patients who are well or poorly perfused in hospitals, hospital-type facilities, mobile, and home environments.

### Comparison to Predicate Devices

The predicate devices used in this filing are:

## 510(k) SUMMARY

- LNCS Oximetry Sensors-Sterile, 510(k) No. K083622
- Reprocessed LNCS Oximetry Sensors, 510(k) No. K083719

The sensors in this filing are substantially equivalent to the predicate sensors in the indication for use/intended use. In addition, they have the same performance and principle of operations, and are subjected to the same processes for sterile (ethylene oxide sterilization) and reprocessing as the respective predicate sensors in the cleared K083622 and K083719 filings. The main difference is that the sensors in this filing have a different interfacing connector in comparison to the predicates.

### Predicate Devices

**Table 1 – Oximetry Sensors, Sterile**

M-LNCS Oximetry Sensors (Single-Use and Sterile)	Predicate LNCS Oximetry Sensors (Single-Use and Sterile) in K083622
M-LNCS Aadx-3: Adult Adhesive Sensor	LNCS Aadx-3
M-LNCS Padx-3: Pediatric Adhesive Sensor	LNCS Padx-3
M-LNCS Inf-3: Infant Adhesive Sensor	LNCS Inf-3
M-LNCS Neo-3: Neonatal Adhesive Sensor	LNCS Neo-3
M-LNCS NeoPt-3: Neonatal Adhesive Sensor	LNCS NeoPt-3
M-LNCS NeoPt-500: Neonatal Adhesive Sensor	LNCS NeoPt-500

**Table 2 – Oximetry Sensors, Reprocessed**

Reprocessed M-LNCS (Single-Use) Sensors	Predicate Reprocessed LNCS (Single-Use) Sensors in K083719
Reprocessed M-LNCS Aadx: Adult Adhesive Sensor	Reprocessed LNCS Aadx
Reprocessed M-LNCS Aadx-3: Adult Adhesive Sensor	Reprocessed LNCS Aadx-3
Reprocessed M-LNCS Padx: Pediatric Adhesive Sensor	Reprocessed LNCS Padx
Reprocessed M-LNCS Padx-3: Pediatric Adhesive Sensor	Reprocessed LNCS Padx-3
Reprocessed M-LNCS Inf: Infant Adhesive Sensor	Reprocessed LNCS Inf
Reprocessed M-LNCS Inf-3: Infant Adhesive Sensor	Reprocessed LNCS Inf-3
Reprocessed M-LNCS Neo: Neonatal Adhesive Sensor	Reprocessed LNCS Neo
Reprocessed M-LNCS Neo-3: Neonatal Adhesive Sensor	Reprocessed LNCS Neo-3
Reprocessed M-LNCS NeoPt: Neonatal Adhesive Sensor	Reprocessed LNCS NeoPt
Reprocessed M-LNCS NeoPt-3: Neonatal Adhesive Sensor	Reprocessed LNCS NeoPt-3
Reprocessed M-LNCS NeoPt-500: Neonatal Adhesive Sensor	Reprocessed LNCS NeoPt-500

**Table 3 - Sensor Specifications for Sterile and Reprocessed M-LNCS Sensors**

	Accuracy Range	Accuracy: Adult/ Pediatric/ Infant	Accuracy: Neonatal
<b>Masimo Technology</b>			
SpO <sub>2</sub> , no motion	70-100%	+ 2%	+ 3%
SpO <sub>2</sub> , motion	70-100%	+ 3%	+ 3%
SpO <sub>2</sub> , low perfusion	70-100%	+ 2%	+ 3%
Pulse rate, no motion	25-240 bpm	+ 3 bpm	+ 3 bpm
Pulse rate, motion	25-240 bpm	+ 5 bpm	+ 5 bpm
Pulse rate, low perfusion	25-240 bpm	+ 3 bpm	+ 3 bpm
<b> Nellcor Technology</b>			
SpO <sub>2</sub> , no motion	70-100%	+ 2%	+ 3%
Pulse rate, no motion	25-240 bpm	± 3 bpm	± 3 bpm

## **510(k) SUMMARY**

### **Test Summary**

The following non-clinical testing was conducted and verifies that the sterile and reprocessed M-LNCS Oximetry Sensors met all design specifications: biocompatibility testing, in-house and laboratory validation testing on the sensors after they have been subjected to ethylene oxide sterilization (EO) or the sensors have been reprocessed and subjected to EO sterilization, performance testing including bench accuracy testing and visual and validated functional testing.

### **Conclusion**

The results of the performance data demonstrate that the sterile M-LNCS and Reprocessed M-LNCS Oximetry Sensors are as safe and effective as the legally marketed predicate devices.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room --WO66-G609  
Silver Spring, MD 20993-0002

Ms. Shelly Harris  
Manager of Regulatory Affairs  
Masimo Corporation  
40 Parker  
Irvine, California 92618

JUL 28 2011

Re: K110723  
Trade/Device Name: Masimo M-LNCS Oximetry Sensors Masimo Reprocessed  
M-LNCS Oximetry Sensors  
Regulation Number: 21 CFR 870.2700  
Regulation Name: Oximeter  
Regulatory Class: II  
Product Code: NLF, DQA, DSA  
Dated: June 14, 2011  
Received: June 28, 2011

Dear Ms. Harris:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Page 2 – Ms. Harris

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours



Anthony D. Watson, B.S., M.S., M.B.A.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K110723

Device Name: Masimo M-LNCS Oximetry Sensors  
Masimo Reprocessed M-LNCS Oximetry Sensors

### Indications For Use:

The M-LNCS Sensors are indicated for the continuous non-invasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO<sub>2</sub>) and pulse rate (measured by an SpO<sub>2</sub> sensor) for use with adult, pediatric, infant, and neonatal patients during both no motion and motion conditions, and for patients who are well or poorly perfused in hospitals, hospital-type facilities, mobile, and home environments.

The Reprocessed M-LNCS Sensors are indicated for the continuous non-invasive monitoring of functional oxygen saturation of arterial hemoglobin (SpO<sub>2</sub>) and pulse rate (measured by an SpO<sub>2</sub> sensor) for use with adult, pediatric, infant, and neonatal patients during both no motion and motion conditions, and for patients who are well or poorly perfused in hospitals, hospital-type facilities, mobile, and home environments.

Prescription Use   X    
(Per 21 CFR 801.109 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(Per 21 CFR 801.109 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)



Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)

Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K110723