

**510(K) SUMMARY**

AUG - 4 2011

**Applicant** IQuum, Inc.  
700 Nickerson Road  
Marlborough, MA 01752  
Tel: 508-970-0099  
Fax: 508-970-0119

**Contact** Lingjun Chen  
Title: Vice President, Operations & Business Development  
Tel: 508-970-0099 ext. 116  
Email: [lingjun@iquum.com](mailto:lingjun@iquum.com)

**Summary Date** May 16, 2011

**Trade Name** Liat™ Influenza A/B Assay, Liat™ Analyzer  
**Common Name** Influenza A, B Panel  
**Regulation Number** 21 CFR 866.3980  
**Classification Name** Respiratory Viral Panel Multiplex Nucleic Acid Assay

**Predicate Device:** Cepheid Xpert® Flu Assay (K103766)

**Intended Use**

The IQuum Liat™ Influenza A/B Assay performed on the Liat™ Analyzer is an automated multiplex real-time RT-PCR assay for the rapid *in vitro* qualitative detection and discrimination of influenza A virus and influenza B virus RNA in nasopharyngeal swab specimens from patients with signs and symptoms of respiratory infection in conjunction with clinical and epidemiological risk factors. The test is intended for use as an aid in the differential diagnosis of influenza A and influenza B in humans and is not intended to detect influenza C.

Negative results do not preclude influenza virus infection and should not be used as the sole basis for treatment or other patient management decisions. Conversely, positive results do not rule-out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.

Performance characteristics for influenza A were established when influenza A/H1 and A/H3 were the predominant influenza A viruses in circulation. When other influenza A viruses are emerging, performance characteristics may vary.

If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent Influenza viruses and sent to state or local health department for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.

## Device Description

The Liat™ Influenza A/B Assay is a rapid, automated in vitro diagnostic test for the detection and differentiation of Influenza type A and type B viral RNA in nasopharyngeal swab (NPS) specimens in universal transport media (UTM) from patients signs and symptoms of respiratory infection. The assay targets a well-conserved region of the matrix gene of Influenza A viral RNA (Inf A target) and non-structural protein (NS) gene of Influenza B (Inf B target). An Internal Process Control (IPC) is also included. The IPC is present to control for adequate processing of the target viruses and to monitor the presence of inhibitors in the RT-PCR reactions.

The Liat Influenza A/B Assay is performed on the lab-in-a-tube technology platform. The system consists of a disposable Liat Influenza A/B Assay Tube and the Liat™ Analyzer. The Liat™ Tube uses a flexible tube as a sample vessel. It contains all required unit dose reagents pre-packed in tube segments, separated by peelable seals, in the order of reagent use. Manipulating a Liat Tube, the Liat Analyzer performs all assay steps from raw sample and report assay result automatically. During the testing process, multiple sample processing actuators of the analyzer compress the Liat Tube to selectively release reagents from tube segments, move the sample from one segment to another, and control reaction volume, temperature and time to conduct sample preparation, nucleic acid extraction, target enrichment, inhibitor removal, nucleic acid elution and real-time RT-PCR. An embedded microprocessor controls and coordinates the actions of these sample processors to perform all required assay processes within the closed Liat Tube. Turnaround time from sample input to result is ~20 minutes.

## Test Operation

A nasopharyngeal swab can be collected following the user institution's standard procedures. For nasopharyngeal swab samples suspended in UTM, the user transfers 100 µl of the UTM sample into Liat Influenza A/B Assay tube. The user then scans the tube barcode to identify the test and scans sample barcode to code the sample ID using the Liat system. The Liat Tube is then inserted into the Liat Analyzer. The analyzer performs all test steps and outputs interpreted results in 20 minutes. A report of the interpreted results can be viewed in the View Results window, and printed directly through a USB connected printer.

No reagent preparation or addition steps are required, other than adding the sample to the Liat tube. Because all the reagents are contained within the Liat assay tube and no sample or reagent needs to be removed from the tube, cross-contamination between samples is minimized.

**Predicate Device Comparison**

Predicate Device: Cepheid Xpert® Flu Assay (K103766)

<b>Item Name</b>	<b>Device: Liat Influenza A/B</b>	<b>Predicate: Cepheid Xpert Flu</b>
Intended Use	<p>The IQuum Liat™ Influenza A/B Assay performed on the Liat™ Analyzer is an automated multiplex real-time RT-PCR assay for the rapid <i>in vitro</i> qualitative detection and discrimination of influenza A virus and influenza B virus RNA in nasopharyngeal swab specimens from patients with signs and symptoms of respiratory infection in conjunction with clinical and epidemiological risk factors. The test is intended for use as an aid in the differential diagnosis of influenza A and influenza B in humans and is not intended to detect influenza C.</p> <p>Negative results do not preclude influenza virus infection and should not be used as the sole basis for treatment or other patient management decisions. Conversely, positive results do not rule-out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.</p> <p>Performance characteristics for influenza A were established when influenza A/H1 and A/H3 were the predominant influenza A viruses in circulation. When other influenza A viruses are emerging, performance characteristics may vary.</p> <p>If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent Influenza viruses and sent to state or local health department for testing. Viral culture should not be attempted in these</p>	<p>The Cepheid® Xpert Flu Assay is an automated, multiplex real-time RT-PCR assay intended for the <i>in vitro</i> qualitative detection and differentiation of influenza A, influenza B and 2009 H1N1 influenza viral RNA. The Xpert Flu Assay uses nasal aspirates/washes and nasopharyngeal swab specimens collected from patients with signs and symptoms of respiratory infection in conjunction with clinical and epidemiological risk factors. The Xpert Flu Assay is intended as an aid in the diagnosis of influenza.</p> <p>Negative results do not preclude influenza virus infection and should not be used as the sole basis for treatment or other patient management decisions.</p> <p>Performance characteristics for influenza A were established during the 2009-2010 influenza season when 2009 H1N1 influenza was the predominant influenza A virus in circulation. When other influenza A viruses are emerging, performance characteristics may vary.</p> <p>If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to state or local health department for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.</p>

	cases unless a BSL 3+ facility is available to receive and culture specimens.	
<b>Item Name</b>	<b>Device: Liat Influenza A/B</b>	<b>Predicate: Cepheid Xpert Flu</b>
Regulation	21 CFR 866.3980	21 CFR 866.3332
Product Code	OCC, OOI	OQW, OCC, OOI
Assay Target	Influenza A Influenza B	Influenza A Influenza B Influenza A subtype 2009 H1N1
Sample Type	Nasopharyngeal Swab	Nasopharyngeal Swab Nasal aspirates/washes
Internal Control	Yes	Yes
Influenza A Viral Target	Matrix gene	Matrix gene
Influenza B Viral Target	Non-structural protein (NSP) gene	Hemagglutinin gene
Assay Instrument	Liat™ Analyzer	GeneXpert instruments
Self-contained System	Integrated PC, software, and touch-screen display	External PC computer and software for running tests and viewing the results
All Assay Reagents Contained in Disposable	No manual reagent addition required	Requires manual dispensing of Binding Reagent into the chamber with the small opening of the Xpert Flu Assay Cartridge
Automated Assay	Yes, sample preparation, amplification, detection and result interpretation	Yes, sample preparation, amplification, detection and result interpretation
Extraction Method	Sample preparation integrated in Liat Tube and Liat Analyzer	Sample preparation integrated in GeneXpert Cartridge and GeneXpert Instrument
Assay Method	RT-PCR for detecting the presence / absence of viral RNA in clinical specimens	RT-PCR for detecting the presence / absence of viral RNA in clinical specimens
Detection Technique	Multiplex assay using different reporter dyes for each target	Multiplex assay using different reporter dyes for each target
Result Interpretation	Automated	Automated
Assay Result	Qualitative	Qualitative
Laboratory Users	CLIA moderate complexity and high complexity laboratory	CLIA moderate complexity and high complexity laboratory
Time-to-result	~20 minutes	~75 minutes

## Substantial Equivalence

### Limit of Detection

The Limit of Detection (LOD) of the Liat Influenza A/B Assay was tested using 3 strains of Influenza A (A/Brisbane/10/2007, A/Brisbane/59/2007, and A/NY/01/2009) and 2 strains of Influenza B (B/Malaysia/2506/04 and B/Florida/04/06). The LOD was determined by limiting dilution studies using these titered viruses. The viruses were spiked into nasopharyngeal swab (NPS) sample matrix, and then tested using the Liat Influenza A/B Assay. The LOD was determined as the lowest log virus concentration that was detected  $\geq 95\%$  of the time (i.e. log concentration at which at least 19 out of 20 replicates tested positive). The LOD for 3 strains of Influenza A were  $10^{-2}$  to  $10^{-1}$  TCID<sub>50</sub>/mL, while those for the 2 strains of Influenza B were  $10^{-3}$  to  $10^{-1}$  TCID<sub>50</sub>/mL.

Virus Strain	LOD (TCID <sub>50</sub> /mL)
A/Brisbane/10/2007	$10^{-1}$
A/Brisbane/59/2007	$10^{-2}$
A/NY/01/2009	$10^{-1}$
B/Florida/04/06	$10^{-1}$
B/Malaysia/2506/04	$10^{-3}$

### Reactivity

Reactivity study evaluates the ability of the Liat Influenza A/B Assay to detect influenza strains representing temporal and geographical diversity. The Liat Influenza A/B Assay was evaluated for reactivity with 22 Influenza A strains and 10 Influenza B strains. Influenza A strains included 8 seasonal Influenza A/H1 strains, 8 seasonal Influenza A/H3 strains, 3 Influenza A 2009 H1N1 strains, 3 swine origin Influenza A strains. The Liat Influenza A/B Assay detected all strains tested.

Influenza Strain	Type / Subtype	Test Concentration	Inf A Result	Inf B Result
A/Brisbane/59/2007	Influenza A, seasonal H1	$8.0 \times 10^{-3}$ TCID <sub>50</sub> /mL	+	-
A/New Caledonia/20/99	Influenza A, seasonal H1	$1.0 \times 10^2$ TCID <sub>50</sub> /mL	+	-
A/Solomon Island/3/2006	Influenza A, seasonal H1	$5.0 \times 10^{-2}$ TCID <sub>50</sub> /mL	+	-
A/Mal/302/54	Influenza A, seasonal H1	$1.0 \times 10^3$ TCID <sub>50</sub> /mL	+	-
A/Denver/1/57	Influenza A, seasonal H1	$5.0 \times 10^2$ TCID <sub>50</sub> /mL	+	-
A/FM/1/47	Influenza A, seasonal H1	$1.0 \times 10^2$ TCID <sub>50</sub> /mL	+	-
A/PR/8/34	Influenza A, seasonal H1	$2.5 \times 10^1$ TCID <sub>50</sub> /mL	+	-
A/Weiss/43	Influenza A, seasonal H1	$2.5 \times 10^3$ TCID <sub>50</sub> /mL	+	-
A/Brisbane/10/2007	Influenza A, seasonal H3	$1.0 \times 10^{-1}$ TCID <sub>50</sub> /mL	+	-

Influenza Strain	Type / Subtype	Test Concentration	Inf A Result	Inf B Result
A/Alice	Influenza A, seasonal H3	$5.0 \times 10^1$ TCID <sub>50</sub> /mL	+	-
A/MRC2	Influenza A, seasonal H3	$1.0 \times 10^2$ TCID <sub>50</sub> /mL	+	-
A/Hong Kong/8/68	Influenza A, seasonal H3	$2.0 \times 10^1$ TCID <sub>50</sub> /mL	+	-
A/Victoria/3/75	Influenza A, seasonal H3	$2.5 \times 10^1$ TCID <sub>50</sub> /mL	+	-
A/Wisconsin/67/05	Influenza A, seasonal H3	$5.0 \times 10^{-1}$ TCID <sub>50</sub> /mL	+	-
A/Port Chalmers/1/73	Influenza A, seasonal H3	$5.0 \times 10^2$ TCID <sub>50</sub> /mL	+	-
A/Aichi/2/68	Influenza A, seasonal H3	$2.0 \times 10^2$ CEID <sub>50</sub> /mL	+	-
A/NY/01/2009	Influenza A, 2009 H1N1	$1.0 \times 10^{-1}$ TCID <sub>50</sub> /mL	+	-
A/NY/02/2009	Influenza A, 2009 H1N1	$2.5 \times 10^{-2}$ TCID <sub>50</sub> /mL	+	-
A/NY/03/2009	Influenza A, 2009 H1N1	$4.0 \times 10^{-1}$ TCID <sub>50</sub> /mL	+	-
A/New Jersey/8/76	Influenza A, H1N1 non 2009	$1.0 \times 10^1$ TCID <sub>50</sub> /mL	+	-
A/Swine/1976/31	Influenza A, H1N1 non 2009	$2.0 \times 10^1$ TCID <sub>50</sub> /mL	+	-
A/Swine/Iowa/15/30	Influenza A, H1N1 non 2009	$2.0 \times 10^2$ TCID <sub>50</sub> /mL	+	-
B/Florida/04/06	Influenza B	$8.0 \times 10^{-2}$ TCID <sub>50</sub> /mL	-	+
B/Malaysia/2506/04	Influenza B	$2.0 \times 10^{-3}$ TCID <sub>50</sub> /mL	-	+
B/Florida/7/04	Influenza B	$5.0 \times 10^{-2}$ TCID <sub>50</sub> /mL	-	+
B/Allen/45	Influenza B	$5.0 \times 10^{-1}$ CEID <sub>50</sub> /mL	-	+
B/GL/1739/54	Influenza B	$2.0 \times 10^1$ TCID <sub>50</sub> /mL	-	+
B/Taiwan/2/62	Influenza B	$5.0 \times 10^{-2}$ TCID <sub>50</sub> /mL	-	+
B/Maryland/1/59	Influenza B	$5.0 \times 10^{-3}$ TCID <sub>50</sub> /mL	-	+
B/Mass/3/66	Influenza B	$1.0 \times 10^1$ TCID <sub>50</sub> /mL	-	+
B/HongKong/5/72	Influenza B	$2.5 \times 10^{-1}$ TCID <sub>50</sub> /mL	-	+
B/Lee/40	Influenza B	$1.0 \times 10^0$ TCID <sub>50</sub> /mL	-	+

Cross Reactivity

Cross reactivity study evaluates the Liat Influenza A/B Assay's potential cross-reactivity with non-influenza respiratory pathogens and other microorganisms with which the majority of the population may have been infected. The Liat assay was evaluated against a panel of 31 human pathogens. Bacteria were tested at  $10^5$ – $10^6$  CFU/mL. Viruses were tested at  $10^3$ – $10^5$  TCID<sub>50</sub>/mL. The Liat Influenza A/B Assay showed no cross reactivity for the tested organisms.

Viral/Bacterial pathogen	Test concentration	Inf A Result	Inf B Result
Adenovirus Type 1	$8.9 \times 10^5$ TCID <sub>50</sub> /mL	–	–
Adenovirus Type 7	$4.5 \times 10^4$ TCID <sub>50</sub> /mL	–	–
Human Coronavirus 229E	$1.4 \times 10^3$ TCID <sub>50</sub> /mL	–	–
Human Coronavirus OC43	$7.9 \times 10^4$ TCID <sub>50</sub> /mL	–	–
Enterovirus	$1 \times 10^5$ TCID <sub>50</sub> /mL	–	–
Human Parainfluenza Type 1	$2.8 \times 10^3$ TCID <sub>50</sub> /mL	–	–
Human Parainfluenza Type 2	$1.4 \times 10^5$ TCID <sub>50</sub> /mL	–	–
Human Parainfluenza Type 3	$1.6 \times 10^5$ TCID <sub>50</sub> /mL	–	–
Measles	$7.9 \times 10^4$ TCID <sub>50</sub> /mL	–	–
Human Metapneumovirus	$7 \times 10^3$ TCID <sub>50</sub> /mL	–	–
Mumps virus	$7.9 \times 10^4$ TCID <sub>50</sub> /mL	–	–
Respiratory syncytial virus type B	$1.4 \times 10^4$ TCID <sub>50</sub> /mL	–	–
Rhinovirus Type 1A	$1.6 \times 10^5$ TCID <sub>50</sub> /mL	–	–
Cytomegalovirus	$4.5 \times 10^4$ TCID <sub>50</sub> /mL	–	–
Epstein Barr virus	$1.9 \times 10^5$ copies/mL	–	–
<i>Bordetella pertussis</i>	$1.8 \times 10^5$ CFU/mL	–	–
<i>Chlamydia pneumoniae</i>	$8 \times 10^4$ TCID <sub>50</sub> /mL	–	–
<i>Corynebacterium sp.</i>	$5.0 \times 10^6$ CFU/mL	–	–
<i>Escherichia coli</i>	$6.6 \times 10^6$ CFU/mL	–	–
<i>Haemophilus influenzae</i>	$3 \times 10^6$ CFU/mL	–	–
<i>Lactobacillus sp.</i>	$1.6 \times 10^6$ CFU/mL	–	–
<i>Legionella pneumophila</i>	$7 \times 10^6$ CFU/mL	–	–
<i>Moraxella catarrhalis</i>	$5.8 \times 10^6$ CFU/mL	–	–
<i>Neisseria meningitidis</i>	$3.2 \times 10^6$ CFU/mL	–	–
<i>Neisseria sp.</i>	$1.8 \times 10^6$ CFU/mL	–	–
<i>Pseudomonas aeruginosa</i>	$1.6 \times 10^6$ CFU/mL	–	–
<i>Staphylococcus aureus</i>	$4.5 \times 10^6$ CFU/mL	–	–
<i>Staphylococcus epidermidis</i>	$6 \times 10^6$ CFU/mL	–	–
<i>Streptococcus pneumoniae</i>	$1.9 \times 10^6$ CFU/mL	–	–
<i>Streptococcus pyogenes</i>	$3.7 \times 10^6$ CFU/mL	–	–
<i>Streptococcus salivarius</i>	$4.3 \times 10^6$ CFU/mL	–	–

Inhibition by other Microorganisms

Interfering microorganism study evaluates whether non-influenza respiratory pathogens and other microorganisms with which the majority of the population may have been infected can interfere in the detection of Influenza A or B by the Liat assay. The panel of 31 human pathogens tested in the cross-reactivity study was tested for potential interference. Bacteria were tested at  $10^5$ – $10^6$  CFU/mL and viruses were tested at  $10^3$ – $10^5$  TCID<sub>50</sub>/mL in the presence of either A/Brisbane/59/2007 or B/Malaysia/2506/04 at 3x LOD concentration in negative NPS matrix. Results show that the presence of the tested microorganisms did not interfere with the detection of Inf A or Inf B.

Pathogen	Pathogen Concentration	A/Brisbane/59/07		B/Malaysia/2506/04	
		Inf A Result	Inf B Result	Inf A Result	Inf B Result
Adenovirus Type 1	8.9×10 <sup>5</sup> TCID <sub>50</sub> /mL	+	–	–	+
Adenovirus Type 7	4.5×10 <sup>4</sup> TCID <sub>50</sub> /mL	+	–	–	+
Human Coronavirus 229E	1.4×10 <sup>3</sup> TCID <sub>50</sub> /mL	+	–	–	+
Human Coronavirus OC43	7.9×10 <sup>4</sup> TCID <sub>50</sub> /mL	+	–	–	+
Enterovirus	1×10 <sup>5</sup> TCID <sub>50</sub> /mL	+	–	–	+
Human Parainfluenza Type 1	2.8×10 <sup>3</sup> TCID <sub>50</sub> /mL	+	–	–	+
Human Parainfluenza Type 2	1.4×10 <sup>5</sup> TCID <sub>50</sub> /mL	+	–	–	+
Human Parainfluenza Type 3	1.6×10 <sup>5</sup> TCID <sub>50</sub> /mL	+	–	–	+
Measles	7.9×10 <sup>4</sup> TCID <sub>50</sub> /mL	+	–	–	+
Human Metapneumovirus	7×10 <sup>3</sup> TCID <sub>50</sub> /mL	+	–	–	+
Mumps virus	7.9×10 <sup>4</sup> TCID <sub>50</sub> /mL	+	–	–	+
Respiratory syncytial virus type B	1.4×10 <sup>4</sup> TCID <sub>50</sub> /mL	+	–	–	+
Rhinovirus Type 1A	1.6×10 <sup>5</sup> TCID <sub>50</sub> /mL	+	–	–	+
Cytomegalovirus	4.5×10 <sup>4</sup> TCID <sub>50</sub> /mL	+	–	–	+
Epstein Barr virus	1.9×10 <sup>5</sup> copies/mL	+	–	–	+
<i>Bordetella pertussis</i>	1.8×10 <sup>5</sup> CFU/mL	+	–	–	+
<i>Chlamydia pneumoniae</i>	8×10 <sup>4</sup> TCID <sub>50</sub> /mL	+	–	–	+
<i>Corynebacterium sp.</i>	5.0×10 <sup>6</sup> CFU/mL	+	–	–	+
<i>Escherichia coli</i>	6.6×10 <sup>6</sup> CFU/mL	+	–	–	+
<i>Haemophilus influenzae</i>	3×10 <sup>6</sup> CFU/mL	+	–	–	+
<i>Lactobacillus sp.</i>	1.6×10 <sup>6</sup> CFU/mL	+	–	–	+
<i>Legionella pneumophila</i>	7×10 <sup>6</sup> CFU/mL	+	–	–	+
<i>Moraxella catarrhalis</i>	5.8×10 <sup>6</sup> CFU/mL	+	–	–	+
<i>Neisseria meningitidis</i>	3.2×10 <sup>6</sup> CFU/mL	+	–	–	+
<i>Neisseria sp.</i>	1.8×10 <sup>6</sup> CFU/mL	+	–	–	+
<i>Pseudomonas aeruginosa</i>	1.6×10 <sup>6</sup> CFU/mL	+	–	–	+
<i>Staphylococcus aureus</i>	4.5×10 <sup>6</sup> CFU/mL	+	–	–	+
<i>Staphylococcus epidermidis</i>	6×10 <sup>6</sup> CFU/mL	+	–	–	+
<i>Streptococcus pneumoniae</i>	1.9×10 <sup>6</sup> CFU/mL	+	–	–	+
<i>Streptococcus pyogenes</i>	3.7×10 <sup>6</sup> CFU/mL	+	–	–	+
<i>Streptococcus salivarius</i>	4.3×10 <sup>6</sup> CFU/mL	+	–	–	+



### Interfering Substances

The Liat Influenza A/B Assay was evaluated with potentially interfering substances that may be encountered in respiratory specimens. Medically and/or physiologically relevant concentrations of potential interferents were tested with 2 influenza A strains and 2 influenza B strains at 3x LOD ( $10^{-1}$ – $10^{-2}$  TCID<sub>50</sub>/mL). Results showed that substances tested did not interfere in the detection of influenza A and B strains.

Potential Interferent	Active Ingredient	Concentration
Mucin Bovine submaxillary gland, type I-S	Purified mucin protein	0.1 mg/ml and 25 mg/ml
Blood	-	5% (v/v)
Nasal spray – Afrin	Oxymetazoline	5% (v/v)
Nasal corticosteroids – Veramyst	Fluticasone	5% (v/v)
Nasal gel – Zicam	Galphimia glauca, Histaminum hydrochloricum, Luffa operculata, Sulphur	5% (v/v)
Throat lozenges, oral anesthetic and analgesic – Cepacol	Benzocaine, Menthol	5 mg/ml
Antibiotic, nasal ointment – Bactroban	Mupirocin	5 mg/ml
Antiviral drug – Relenza	Zanamivir	5 mg/ml
Antiviral drug – Tamiflu	Oseltamivir	7.5 mg/ml
Antimicrobial, systemic	Tobramycin	4 µg/ml

### Inter-lot Precision

Inter-lot precision assesses the Liat tube manufacturing process as a potential source of assay variability. Three lots of Liat Influenza Assay tubes were tested with 1 Influenza A strain at C95 (LOD, n=180), C100 (4xLOD, n=60) and C0 (negative NPS matrix, n=60). Results from all runs agreed with expected results. Imprecision in Ct from inter-lot variability was <1.8%.

Inf A	n	% Agreement	Avg Ct	Stdev Inter-lot	%CV Inter-lot
C100	60	100%	29.6	0.49	1.6%
C95	180	100%	32.0	0.59	1.8%
C0	60	100%	42.0	0.00	0.0%

### Reproducibility

Reproducibility study assesses the total variability of the Liat Influenza A/B Assay across operators, study sites, testing days, Liat Analyzers, and Liat assay tube lots. The Liat assay was evaluated at 3 sites, including 2 CLIA waived sites. Two operators at each of the 3 sites tested an 8 member reproducibility panel in triplicate on 5 different days, for a total of 720 runs (8 panel members × 3 replicates × 2 operators × 5 days × 3 sites). Fifteen (15) Liat Analyzers and 3 Liat Influenza A/B Assay tube lots were used. The reproducibility panel comprises a high negative, a low positive and a medium positive of each of Influenza A and B, and the assay positive and negative controls. Percent agreement with expected result, mean Ct, and Ct %CV for each site are show in tables below. %CV for Influenza A ranged between 1.3% and 3.8% and that for Influenza B ranged between 1.1% and 3.4%. Total percent agreement was ≥99.9%.

**Influenza A Reproducibility**

Sample	Site 1			Site 2			Site 3			Total	
	Agreement w/ expected result	Avg Ct	Total %CV	Agreement w/ expected result	Avg Ct	Total %CV	Agreement w/ expected result	Avg Ct	Total %CV	Agreement w/ expected result	95% CI
Positive Control	30/30	32.0	1.4%	30/30	32.0	2.2%	30/30	31.5	1.3%	90/90 (100%)	95.9%-100.0%
Negative Control	30/30	39.0	0.0%	30/30	39.0	0.0%	30/30	39.0	0.0%	90/90 (100%)	95.9%-100.0%
Flu A High Negative	30/30	39.0	0.0%	30/30	39.0	0.0%	30/30	39.0	0.0%	90/90 (100%)	95.9%-100.0%
Flu A Low Positive	30/30	32.1	3.8%	30/30	32.0	2.1%	30/30	31.8	1.8%	90/90 (100%)	95.9%-100.0%
Flu A Med. Positive	30/30	29.9	2.0%	30/30	29.6	2.0%	30/30	29.3	1.9%	90/90 (100%)	95.9%-100.0%
Flu B High Negative	30/30	39.0	0.0%	30/30	39.0	0.0%	30/30	39.0	0.0%	90/90 (100%)	95.9%-100.0%
Flu B Low Positive	30/30	39.0	0.0%	30/30	39.0	0.0%	30/30	39.0	0.0%	90/90 (100%)	95.9%-100.0%
Flu B Med. Positive	30/30	39.0	0.0%	30/30	39.0	0.0%	30/30	39.0	0.0%	90/90 (100%)	95.9%-100.0%
Total Agreement	240/240 (100%)			240/240 (100%)			240/240 (100%)			720/720 (100%)	99.5%-100.0%

**Influenza B Reproducibility**

Sample	Site 1			Site 2			Site 3			Total	
	Agreement w/ expected result	Avg Ct	Total %CV	Agreement w/ expected result	Avg Ct	Total %CV	Agreement w/ expected result	Avg Ct	Total %CV	Agreement w/ expected result	95% CI
Positive Control	30/30	31.2	1.5%	30/30	31.0	2.2%	30/30	30.6	1.1%	90/90 (100%)	95.9%-100.0%
Negative Control	30/30	39.0	0.0%	30/30	39.0	0.0%	30/30	39.0	0.0%	90/90 (100%)	95.9%-100.0%
Flu A High Negative	30/30	39.0	0.0%	30/30	39.0	0.0%	30/30	39.0	0.0%	90/90 (100%)	95.9%-100.0%
Flu A Low Positive	30/30	39.0	0.0%	30/30	39.0	0.0%	30/30	39.0	0.0%	90/90 (100%)	95.9%-100.0%
Flu A Med. Positive	30/30	39.0	0.0%	30/30	39.0	0.0%	30/30	39.0	0.0%	90/90 (100%)	95.9%-100.0%
Flu B High Negative	30/30	39.0	0.0%	29/30	38.9	1.7%	30/30	39.0	0.0%	89/90 (98.9%)	94.0%-99.8%
Flu B Low Positive	30/30	31.3	2.3%	30/30	31.0	3.4%	30/30	30.9	1.7%	90/90 (100%)	95.9%-100.0%
Flu B Med. Positive	30/30	29.4	2.6%	30/30	29.5	3.0%	30/30	28.6	3.3%	90/90 (100%)	95.9%-100.0%
Total Agreement	240/240 (100%)			239/240 (99.6%)			240/240 (100%)			719/720 (99.9%)	99.2%-100.0%

### Performance using Fresh vs. Frozen Samples

The Liat Influenza A/B Assay was tested by comparing its performance using fresh and frozen specimens. One influenza A strain and one influenza B strain were individually spiked into NPS matrix at different viral loads, including levels near LOD and levels reflecting the clinical range. For each strain, 60 samples were tested immediately while another 60 samples were frozen at -80°C for 7 days, thawed and then tested. Fresh and frozen samples demonstrated 100% detection across all levels of viral load, demonstrating that the Liat Influenza A/B Assay had equivalent performance for fresh and frozen samples.

### Clinical Agreement

Three sites participated in the clinical study for the Liat Influenza A/B Assay. A total of 615 clinical specimens were tested including 435 prospectively collected samples and 180 retrospective samples. Prospective samples were collected from patients with signs and symptoms of influenza in the Eastern and Southwestern US from 12 February 2009 to 26 March 2009. Viral culture and IFA staining was used as the reference method for these prospectively collected samples. Discordance results were investigated using PCR and bi-directional sequencing.

Due to the low prevalence of influenza during the collection period as well as the subsequent emergence of the 2009 H1N1 virus, retrospective samples collected between 2008-2010, including 2009 H1N1 samples from the 2009-2010 flu season, were tested. Reference testing for these samples was performed by PCR and bi-directional sequencing based on published methods (Ghedini et al. 2005; Ghedini et al. 2009; World Health Organization Sequencing Primers and Protocols, 2009). Of the 180 samples, 1 sample was indeterminate by PCR/sequencing, and remained indeterminate upon retest; this sample was excluded from the analysis:

Clinical sample testing was performed at 3 sites in April-May 2011 using 30 Liat Analyzers in total. Tables below summarize the clinical performance of the Liat Influenza A/B Assay.

### **Liat Influenza A/B Assay Clinical Performance – Prospective Samples**

<b>Influenza A</b>		Viral Culture				%	95% CI
		Positive	Negative	Total			
Liat	Positive	34	13 <sup>a</sup>	47	Sensitivity	100.0%	(89.8% - 100.0%)
	Negative	0	388	388	Specificity	96.8%	(94.5% - 98.1%)
Total		34	401	435			

<sup>a</sup> Of 13 false positive samples, 8 were Influenza A positive by PCR/sequencing, 4 were negative by PCR/sequencing, and 1 was indeterminate by PCR/sequencing due to low sequence quality score.

<b>Influenza B</b>		Viral Culture				%	95% CI
		Positive	Negative	Total			
Liat	Positive	30	24 <sup>b</sup>	54	Sensitivity	100.0%	(88.6% - 100.0%)
	Negative	0	381	381	Specificity	94.1%	(91.3% - 96.0%)
Total		30	405	435			

<sup>b</sup> Of 24 false positive samples, 13 were Influenza B positive by PCR/sequencing, 3 were negative by PCR/sequencing, and 8 samples were indeterminate by PCR/sequencing due to low sequence

quality score.

### Liat Influenza A/B Assay Clinical Performance – Retrospective Samples

Influenza A		PCR/Sequencing		Total
		Positive	Negative	
Liat	Positive	74 <sup>c</sup>	3	77
	Negative	0	102	102
Total		74	105	179

	%	95% CI
Positive Agreement	100.0%	(95.1% - 100.0%)
Negative Agreement	97.1%	(91.9% - 99.0%)

<sup>c</sup> Of 74 Influenza A positive retrospective samples, 44 were 2009 H1N1 positive, and 20 were A/H3 positive by PCR/sequencing. The 10 remaining samples were verified to be Influenza A however, no further subtyping analysis was done.

Influenza B		PCR/Sequencing		Total
		Positive	Negative	
Liat	Positive	7	1	8
	Negative	0	171	171
Total		7	172	179

	%	95% CI
Positive Agreement	100.0%	(64.6% - 100.0%)
Negative Agreement	99.4%	(96.8% - 99.9%)

### Conclusion

The results of the analytical and clinical performance studies submitted in this premarket notification are complete and demonstrate that the Liat Influenza A/B assay on the Liat Analyzer are substantially equivalent to the predicate device.



Food and Drug Administration  
10903 New Hampshire Avenue  
Silver Spring, MD 20993

IQuum Inc.  
c/o Lingjun Chen  
Vice President, Operations and Business Development  
700 Nickerson Road  
Marlborough, MA 01752

AUG - 4 2011

Re: K111387

Trade/Device Name: Liat Influenza A/B Assay, Liat Analyzer  
Regulation Number: 21 CFR§ 866.3980  
Regulation Name: Respiratory Viral Panel Multiplex Nucleic Acid Assay  
Regulatory Class: Class II  
Product Code: OCC, OOI  
Dated: May 16, 2011  
Received: May 18, 2011

Dear Dr. Chen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into class II (Special Controls), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket

notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Parts 801 and 809), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Sally A. Hojvat, M.Sc., Ph.D.  
Director  
Division of Microbiology Devices  
Office of *In Vitro* Diagnostic Device  
Evaluation and Safety  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K111387

Device Name: Liat™ Influenza A/B Assay, Liat™ Analyzer

### Indications for Use:

The IQuum Liat™ Influenza A/B Assay performed on the Liat™ Analyzer is an automated multiplex real-time RT-PCR assay for the rapid *in vitro* qualitative detection and discrimination of influenza A virus and influenza B virus RNA in nasopharyngeal swab specimens from patients with signs and symptoms of respiratory infection in conjunction with clinical and epidemiological risk factors. The test is intended for use as an aid in the differential diagnosis of influenza A and influenza B in humans and is not intended to detect influenza C.

Negative results do not preclude influenza virus infection and should not be used as the sole basis for treatment or other patient management decisions. Conversely, positive results do not rule-out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.

Performance characteristics for influenza A were established when influenza A/H1 and A/H3 were the predominant influenza A viruses in circulation. When other influenza A viruses are emerging, performance characteristics may vary.

If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent Influenza viruses and sent to state or local health department for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF  
NEEDED)

---

Concurrence of CDRH, Office of *In Vitro* Diagnostics (OIVD)

Tamara Feldblum  
Division Sign-Off

Office of In Vitro Diagnostic Device  
Evaluation and Safety

510(k) K111387