

K111836

Ultradent Products, Inc.  
Premarket Submission for PermaShade™ Veneer Cement  
Traditional 510(k)

SEP 16 2011

## Section 5: 510(k) Summary

This summary of the Traditional 510(k) substantial equivalence information is being submitted in accordance with the requirements of 21 CFR 807.92.

### I. Applicant's Name and Address

Ultradent Products, Inc.  
505 West 10200 South  
South Jordan, UT 84095

Contact Person: Corey Jaseph, RAC  
Title: Regulatory Affairs Manager  
Telephone: 800-552-5512 x4420, 801-553-4420  
FAX: 801-553-4609  
Date Summary Prepared: 15 SEP 2011

### II. Name of the Device

Trade Name: PermaShade Veneer Cement  
Common Name: Dental Cement  
Device Classification: 2  
Classification Product Code: EMA  
Regulation No. 872.3275 (dental cement)

### III. Legally Marketed Predicate Devices to Which Equivalence is Claimed

The predicate devices are PermaFlo, manufactured by Ultradent (K974413), and Variolink Veneer, manufactured by Ivoclar Vivadent (K931309).

### IV. Device Description:

PermaShade Veneer Cement is a light-cured, moderately viscous, permanent resin cement. It is 66.5% filled by mass, is available in four shades, and has very low shrinkage and color change properties.

**V. Statement of intended use:**

PermaShade Veneer Cement can be used for the permanent cementation of porcelain, zirconia, composite, and other indirect veneers.

**VI. Comparison of technological characteristics**

**Table 5-1: Substantial equivalence comparison**

Characteristic	PermaFlo Ultradent K974413	Variolink Veneer Ivoclar Vivadent K931309	PermaShade Veneer Cement Ultradent NEW DEVICE
<b>Intended Use</b>	PermaFlo can be used for: 1. Anterior and posterior restorations, such as Class I, II, III, IV, and V 2. Luting translucent inlays and onlays. 3. Direct veneers and other restorative procedures. 4. Restoring missing subgingival tooth structure prior to endodontic procedures (this is referred to as the "donut" technique)	Adhesive luting of ceramic and composite restorations	Permanent cementation of indirect porcelain, zirconia, and composite veneers
<b>Intended user</b>	Dental professional	Same	Same
<b>Composition of materials</b>	Dimethacrylate-based, filled resin	Same	Same
<b>Type of cure</b>	Photochemical	Same	Same
<b>Percent Fill</b>	68% (m/m)	47 – 66% (m/m)	66% (m/m)
<b>Filler Particle Size</b>	Average 700 nm	40 – 300 nm	13 – 700 nm
<b>Delivery system</b>	Syringe	Same	Same
<b>Multiple Shades</b>	Yes	Same	Same
<b>Physical properties</b>	Tested and characterized according to ISO 4049	Same	Same
<b>Biocompatibility</b>	Tested per ISO 10993-1 and ISO 7405	Same	Same

PermaShade Veneer Cement was designed to be a strong and esthetic option with high color stability for cementing indirect veneers. PermaShade Veneer Cement is a similar material used in the same way by the same types of users as the identified predicate devices, introducing no new safety or efficacy questions. It has been tested against the predicates listed above as outlined in ISO 4049 and demonstrated equivalent in vitro

performance. Biocompatibility testing shows that the product is safe when used as instructed by a dental professional. In summary, this submission demonstrates that PermaShade Veneer Cement is safe and effective and performs equivalently to the identified predicates for its intended use.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room -WO66-G609  
Silver Spring, MD 20993-0002

Ms. Corey Jaseph  
Regulatory Affairs Manager  
Ultradent Products, Incorporated  
505 West 10200 South  
South Jordan, Utah 84095

SEP 16 2011

Re: K111836  
Trade/Device Name: PermaShade Veneer Cement  
Regulation Number: 21 CFR 872.3275  
Regulation Name: Dental Cement  
Regulatory Class: II  
Product Code: EMA  
Dated: June 28, 2011  
Received: June 29, 2011

Dear Ms. Jaseph:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Handwritten signature of Anthony D. Watson, consisting of a stylized 'A' followed by 'D. Watson' and the word 'for'.

Anthony D. Watson, B.S., M.S., M.B.A.  
Director  
Division of Anesthesiology, General Hospital,  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Section 4: Statement of Indications for Use**

510(k) Number (if known): K111836

Device Name: PermaShade Veneer Cement

Indications for Use: PermaShade Veneer Cement can be used for permanent cementation of porcelain, zirconia, composite and other indirect veneers.

Prescription Use X  
(Part 21 CFR 801 Subpart D)

AND/OR Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of 1

(Posted November 13, 2003)

Susan Runner  
(Division Sign-Off)  
Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K111836