

AUG 24 2011

SECTION 7 - 510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

Submitter Name & Address: Corventis, Inc.
1410 Energy Park Drive, Suite 1
St. Paul, MN 55108

Contact Person: Michele Chin-Purcell, Ph.D.
651-925-3803 (phone)
651-389-3251 (fax)

Trade/Proprietary Name: NUVANT® Mobile Cardiac Telemetry System

Common/Usual Name: Mobile Cardiac Telemetry System

Classification Name: Arrhythmia Detector and Alarm
(21 CFR 870.1025, Product Code DSI)
Patient Physiological Monitor (with arrhythmia
detection)
(21 CFR 870.1025, Product Code MHX)

Class: Class II, Special Controls

510(k): Abbreviated 510(k)

Date Prepared: July 5, 2011

Predicate Devices:

1. NUVANT Mobile Cardiac Telemetry System, Corventis, Inc. cleared by FDA under 510(k) number K091971; 21 CFR 870.1025, DSI "Arrhythmia Detector and Alarm", and 21 CFR 870.1025, MHX "Patient Physiological Monitor (with arrhythmia detection)"
2. CardioNet ECG Monitor with Arrhythmia Detection Model CN1005, CardioNet Inc., cleared by FDA under K072558; 21 CFR 870.1025, DSI "Arrhythmia Detector and Alarm"

BACKGROUND INFORMATION

The following is the subject device of this Abbreviated 510(k):

- The NUVANT Mobile Cardiac Telemetry System (hereafter referred to as "subject NUVANT")

This submission proposes to add clarifying language (e.g., listing the type of non-lethal arrhythmias) to the predicate NUVANT indications for use. This clarifying language is similar to that used in the predicate device from CardioNet. This change is for clarification and does not change to the intended therapeutic, diagnostic, prosthetic, or surgical use of the device, and does not affect the safety and effectiveness of the device when used as labeled [807.92(a)(5)].

The intent of the subject devices remains the same as that of the predicate devices, which is to monitor ECG and heart rate.

INDICATION FOR USE STATEMENT

The following are the proposed indications for use in this submission:

The NUVANT Mobile Cardiac Telemetry (MCT) System is intended to continuously measure, record, and periodically transmit physiological data. The System is indicated for those patients who require monitoring for the detection of non-lethal cardiac arrhythmias such as, but not limited to, supraventricular tachycardias (e.g. atrial fibrillation, atrial flutter, paroxysmal SVTs), ventricular ectopy, bradyarrhythmias and conduction disorders.

The NUVANT system model monitors, derives and displays:

- ECG
- Heart rate

TECHNOLOGICAL CHARACTERISTICS AND SUBSTANTIAL EQUIVALENCE

The NUVANT MCT System is wearable, wireless arrhythmia detection system that is used by clinicians to identify suspected cardiac arrhythmias. In combination with interpretation services provided by learned intermediaries in the Corventis Monitoring Center as well as online review of data (for prescribing physicians only), NUVANT MCT enables arrhythmia detection for up to 7.5 days for each PiiX application.

The NUVANT system components are:

- PiiX® (aka: Adherent Device) - a patient-worn device which is applied to the patient's torso. It contains the ECG electrodes for recording ECG and heart rate data.
- Patient Trigger Magnet - used by the patient to manually trigger the ECG collection when he/she experiences symptoms.
- zLink® (aka: Gateway) – hand-held device that receives information from the PiiX and transmits it to the Corventis Server via cellular technology.
- Server - The Server receives sensor data from the PiiX via zLink. ECG and heart rate are presented to learned intermediaries, Corventis cardiographic technicians, who prepare and deliver the information to prescribing physicians

The communication between the PiiX and the zLink is enabled via BlueTooth Technology. Sensor data and ECGs collected by the PiiX are transmitted to the Server via zLink.

The system components of the subject NUVANT are identical to the system components of the predicate NUVANT.

PREDICATE DEVICES

Two (2) predicate devices have been identified:

1. NUVANT Mobile Cardiac Telemetry System, Corventis, Inc. cleared by FDA under 510(k) number K091971; 21 CFR 870.1025, DSI "Arrhythmia Detector and Alarm", and 21 CFR 870.1025, MHX "Patient Physiological Monitor (with arrhythmia detection)"
2. CardioNet ECG Monitor with Arrhythmia Detection Model CN1005, CardioNet Inc., cleared by FDA under K072558; 21 CFR 870.1025, DSI "Arrhythmia Detector and Alarm"

The following table summarizes the compared features and the corresponding predicate devices.

Features comparison and corresponding predicate devices

Features being compared	Subject Devices	Predicate Devices
All technical features	Subject NUVANT	Predicate NUVANT (K091971) CardioNet CN1005 (K072558)
Indications for Use	Subject NUVANT	Predicate NUVANT (K091971) CardioNet CN1005 (K072558)

SUMMARY OF PERFORMANCE TESTING

Both the subject and predicate NUVANT MCT meet the requirements of following performance standards in accordance with *FDA Class II Special Controls Guidance Document: Arrhythmia Detector and Alarm*.

IEC 60601-1

Medical Electrical Equipment – Part 1: General Requirements for Safety, 1988;
Amendment 1, 1991-11, Amendment 2, 1995

IEC 60601-1-2

Medical Electrical Equipment - Part 1-2: General Requirements for Safety -
Collateral standard: Electromagnetic Compatibility - Requirements and Tests
(Edition 2:2001 with Amendment 1:2004; Edition 2.1 (Edition 2:2001
consolidated with Amendment 1:2004))

AAMI/ANSI EC38

Medical electrical equipment - Part 2-47: Particular requirements for the safety,
including essential performance, of ambulatory electrocardiographic systems,
2007.

AAMI/ANSI EC57

Testing and reporting performance results of cardiac rhythm and ST-segment
measurement algorithms, 1998/(R) 2008

CONCLUSION

The modifications proposed in the subject NUVANT System make clarifications to the indications for use and do not change the fundamental scientific technology or use of the devices. As supported by the descriptive information and the bench tests, it is concluded that the NUVANT Mobile Cardiac Telemetry (NUVANT) Systems are as safe and as effective as the predicate devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room -WO66-G609
Silver Spring, MD 20993-0002

Corventis, Inc.
c/o Michele Chin-Purcell, Ph.D.
Senior Director, Quality and Regulatory Affairs
1410 Energy Park Drive, Suite 1
St. Paul, MN 55108

AUG 24 2011

Re: K111917
Trade/Device Name: Nuvant Mobile Cardiac Telemetry
Regulatory Number: 21 CFR 870.1025
Regulation Name: Patient Physiological Monitor (with arrhythmia detection or alarms)
Regulatory Class: II (two)
Product Code: 74 MHX
Dated: July 5, 2011
Received: July 6, 2011

Dear Dr. Chin-Purcell:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Page 2 – Michele Chin-Purcell, Ph.D.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Bram D. Zuckerman, M.D.", is written over a horizontal line. The signature is fluid and cursive.

Director
Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

K111917

Indications for Use

510(k) Number (if known): N/A

Device Name: NUVANT® Mobile Cardiac Telemetry

Indications for Use:

The NUVANT Mobile Cardiac Telemetry (MCT) System is intended to continuously measure, record and periodically transmit physiological data. The System is indicated for those patients who require monitoring for the detection of non-lethal cardiac arrhythmias such as, but not limited to, supraventricular tachycardias (e.g. atrial fibrillation, atrial flutter, paroxysmal SVTs), ventricular ectopy, bradyarrhythmias and conduction disorders.

The NUVANT system model monitors, derives and displays:

- ECG
- Heart Rate

Prescription Use X
(Part 21 CFR 801 Subpart D)


AND/OR

Over-The-Counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of 1


(Division Sign-Off)
Division of Cardiovascular Devices

510(k) Number K111917